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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the signature of Huang Xuehua
of JIANGSU BONSS MEDICAL TECHNOLOGY CO.,LTD. and the
seal of the said company on the annexed DOCUMENT are
genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature: Cheng Jia

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(Date: Oct. 19, 2020)

证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

BONSS

Approved by
Director
Huan Xuehua


OPERATION MANUAL

Surgical radio frequency plasma system ARS

execution versions

ARS600, ARS700, ARS800

manufacturer JIANGSU BONSS MEDICAL TECHNOLOGY CO. LTD, China

Version 01
September, 2019

Content

Table of contents

1. Name of medical device, composition	3
2. Manufacturer	4
3. Purpose	4
4. Performance characteristics	5
5. Risks, contra-indications and side effects	6
6. Technical characterisites	9
7. Installation order	12
8. TYPES OF ELECTRODES	25
9. LIST OF MATERIALS COMING INTO CONTACT WITH AN ORGANISM	26
10. SAFE COMBINATIONS with other products	28
11. SOFTWARE	29
12. WORKING AND ALARM SIGNALS	30
13. Cleaning, disinfection and maintenance	31
14. Warranty obligations	31
15. Conditions of storage and transportation	31
16. MARKING AND PACKING	33
17. OPERATION	41
18. Environmental protection	43
19. ELECTROMAGNETIC COMPATIBILITY (EMC)	43
20. Service life, WARRANTY OBLIGATIONS	47
21. RECLAMATIONS	47
22. Disposal	47
23. INTERNATIONAL STANDARDS	47

ATTENTION

This device complies with the requirements of MDD 93/42/EEC.

The device can generate and use radio frequency energy. Improper installation may cause interference for surrounding equipment. The presence of interference is checked by turning the device on / off.

1. NAME OF MEDICAL DEVICE, COMPOSITION

Surgical radio frequency plasma system ARS

1.1. Composition of medical device. Delivery set

Surgical radio frequency plasma system

execution versions:

I. Surgical radio frequency plasma system ARS600 consisting of:

1. Controller ARS600;
2. Power supply cord;
3. Shift pedals AJ120;
4. Fuse 5A – 2 pcs.,
5. Operation manual;
6. Electrode AC301, 301Q2, (not more than 100 pcs.),
7. Electrode AC302, 302Q4, (not more than 100 pcs.),
8. Electrode MC302, 302Q5, (not more than 100 pcs.),
9. Electrode AC306, 306Q1, (not more than 100 pcs.),
10. Electrode AC306, 306Q2, (not more than 100 pcs.),
11. Electrode AC 404, 404Q1, (not more than 100 pcs.),
12. Electrode AC404, 404Q2, (not more than 100 pcs.),
13. Electrode AC405, 405Q2, (not more than 100 pcs.),
14. Electrode AC405, 405Q3, (not more than 100 pcs.).

II. Surgical radio frequency plasma system ARS700 consisting of

1. Controller ARS700,
2. Power supply cord,
3. Shift pedals AJ120,
4. Stream unit,
5. Cable of stream unit,
6. Fuse SA – 2 pcs.,
7. Operation manual,
8. Electrode AC304, A4845S, (not more than 100 pcs.),
9. Electrode AC305, A4855S, (not more than 100 pcs.),
10. Electrode AC401, A5874, (not more than 100 pcs.),
11. Electrode AC401, A5875, (not more than 100 pcs.),
12. Electrode AC402, A8872, (not more than 100 pcs.),
13. Electrode AC403, A7070, (not more than 100 pcs.),
14. Electrode AC403, A7171, (not more than 100 pcs.).

III. Surgical radio frequency plasma system ARS800 consisting of:

1. Controller ARS800,
2. Power supply cord,
3. Shift pedals AJ120,
4. Fuse 5A – 2 pcs.,
5. Operation manual,
6. Electrode AC307, 106Q1, (not more than 100 pcs.),
7. Electrode AC307, 106Q2, (not more than 100 pcs.).

2. MANUFACTURER

JIANGSU BONSS MEDICAL TECHNOLOGY CO. LTD, Building #7, No.898, China Medical City Avenue, Hailing District, Taizhou City, 225316 Jiangsu P.R. China
 Tel. +86-(0)523-86813258
 Fax: +86-(0)523-86813258
sales@bonss.com.cn

Place of manufacturing

JIANGSU BONSS MEDICAL TECHNOLOGY CO. LTD, Building #7, No.898, China Medical City Avenue, Hailing District, Taizhou City, 225316 Jiangsu P.R. China.
 Tel. +86-(0)523-86813258
 Fax: +86-(0)523-86813258
sales@bonss.com.cn

Authorized representative in the Russian Federation

«SMK-Trade» Limited Liability Company
 «SMK-Trade» LLC
 Russia, 214010, Smolensk region, city of Smolensk, Sedov st., 42
 Tel. +7 905 163 58 60
aleksey.zhdanyuk@yandex.ru

3. PURPOSE

Medical device - Surgical radio frequency plasma system ARS (hereinafter referred to as the system, surgical system, device) is intended for ablation, resection and coagulation of soft tissues and hemostasis of blood vessels with the energy of low-temperature plasma in ENT surgery, orthopedics, spine surgery, urology and gynecology.

Potential consumer

Only a qualified surgeon is allowed to operate the product.
 Before performing any surgical procedures with the system, the surgeon must complete a special training course in working with electrosurgical equipment.
 The user must be experienced in performing electrosurgical procedures.

Scope

- ENT surgery,
- Orthopedic surgery
- Spine surgery
- Gynecology
- Urology

3.1. Indications for use


The surgical system can be used for male and female patients from the age of 3 years.

- The need to remove tonsils, adenoids
- Ablation, coagulation and removal of pathologically altered tissues in ENT, such as the nasal cavity, soft palate, laryngeal, vocal cord, tongue base, nasal skull base
- Intra-articular changes of traumatic and degenerative-dystrophic nature

- Herniation and protrusion of intervertebral discs in the cervical and lumbar spine (nucleoplasty method)
- Pathological changes in female genital organs
- Pathological changes in the urinary system

4. PERFORMANCE CHARACTERISTICS

The system is a radio frequency plasma electrosurgical device with bipolar and multipolar functions. There are two modes of system operation:

 Yellow keys of ablation power on the controller panel	1. Ablation mode «AB-LATE» - for resection and ablation, the power of mode is set using the keys of yellow control panel of the controller and is activated by pressing the yellow foot pedal AJ120.	 Shift pedals AJ120 When the pedal is depressed, a sound signal is heard (different for the ablation mode and coagulation mode. (See the section "Working and alarm signals").
 Blue keys for coagulation power on the controller panel	2. Coagulation mode «CO-AG» - for coagulation (hemostasis of blood vessels), the power of which is set using the buttons on the blue control panel of controller and activated by pressing the blue foot pedal AJ120.	

The LED digital display shows the value of ablation and coagulation power in the format from 1 to 9. The procedure time is displayed in the format from 0 to 9, which means 0-900 ms, with an error of +/- 200 ms.

The controller has a function for recognizing the type of electrode and starts counting the activation time from 100 ms.

A green "Power" LED indicates that the power supply is on. When the electrode is connected, the green indicator on the display panel lights up; if the electrode is not connected, the red indicator lights up.

Operating principle

The surgical system generates 100 kHz radio frequency energy.


The basic principle of operation of all electrosurgical devices is to create a heat generation effect. In electrical surgery, heat is generated when charged particles pass through human tissue that resist flow. The process is called "resistive heating". With a sharp heating, the intracellular fluid medium reaches the boiling process, the formation of steam, which entails an increase in intracellular pressure and destruction of cell wall. Delamination of tissue is observed.

The Surgical Radiofrequency Plasma ARS System is a product of cold plasma surgery (coblation technology). Coblation, unlike traditional radiofrequency surgery, heats tissues by passing low-frequency electric currents through the electrolyte medium, causing a cloud of ionized high-energy

particles (plasma). Plasma ions produce heating of tissues to very low temperatures (40-70°C), but they have energy that can destroy biological bonds.

A low-temperature plasma layer is formed in a limited space of conducting medium (in the intercellular fluid of organ or in a saline solution irrigating the electrode) - between the working head (tip) of electrode and organ tissue.

In the plasma ablation mode, saline (NaCl) is used as the working fluid. As an electric current passes through the saline solution, the ions are accelerated and create a thin layer of plasma energy around the electrodes. The thin layer of plasma is composed of massive charged particles that can generate enough energy to break down organic molecular bonds within tissue and cause tissue to dissolve rapidly at relatively low temperatures.

	The RF energy emitted from the activated electrode forms a precisely focused plasma envelope of conductive saline or tissue fluid ions (which is also an electrolyte). Plasma, consisting of massive charged particles, has enough energy to destroy the organic molecular bonds of tissues. The destruction of tissues to the molecular and atomic level occurs almost instantly at relatively low temperatures.
	When RF energy is applied to tissue (including blood) around the electrode, Joule heat is generated and an electromagnetic wave effect is generated to immediately coagulate (clot) the tissue protein (blood).

When compared with standard burning, which occurs at temperatures between 400 and 600 ° C, the coblation temperature, along with simultaneous irrigation with saline, can provide very little thermal effect on the surrounding tissues. Where the plasma layer meets tissue, there is sufficient energy to break molecular bonds, resulting in molecular dissociation. Since this effect is achieved at a temperature of approximately 40 to 70 ° C, thermal damage to the tissue is minimized.

Coblation creates a relatively bloodless surgical field that provides accurate visualization of surgical site. These factors help to improve accuracy and potentially reduce tissue damage and subsequent fibrosis for the surrounding structures. Coblation involves the passage of bipolar current through saline, which doesn't spark in an oxygen-rich environment. This significantly reduces the risk of airway inflammation compared to unipolar electrocoagulation methods.

Each controller version is used only with the electrodes included in its composition. Some types of electrodes have the function of saline irrigation, working with an electrolyte dispenser. In the ARS700 version, this function is performed by the flow control unit. The flow control unit is used to control irrigation through the electrode tubes. It works in sync with the generator. The unit can be turned on or off automatically when the generator is activated or turned off.

5. RISKS, CONTRAINDICATIONS, SIDE EFFECTS

Risk analysis and management was carried out in accordance with EN ISO 14971: 2012.

Risk sources were considered: energy hazards, biological and chemical hazards, operational hazards (functional and user errors), information hazards.

- All of the above risks have been addressed and mitigated as low as possible. Residual risks are controlled and acceptable. They are indicated in the operation manual in the paragraphs "Warnings" and "Precautions".

- Risks that can't be mitigated by risk management measures, such as prohibiting the use of surgical system for patients with an implanted pacemaker or other electronic implant, are listed in the section "Contraindications".

Contraindications for use

- Patients using a pacemaker or other electronic implants are prohibited from performing surgical operations using the surgical radio frequency plasma system.
- The use of surgical system is contraindicated for procedures in the absence of electrically conductive medium.

Possible side effects

- Post-operative pain
- Bleeding

Warnings

Attention! The system generates radio frequency energy and may interfere with other electrical appliances.

- Read the Operation Manual carefully before using the surgical system. Observe all warnings and precautions.
- Activated or hot (immediately after use) electrodes can cause a fire. It is forbidden to activate the electrodes in contact with the surface of System and its components, as well as with flammable materials: gauze or cloth coatings, as this may lead to a fire.
- Don't touch the working tip (head) of activated electrode.
- After disconnecting the controller from the mains, the electrode tips may remain hot enough to cause burns.
- Don't activate or move electrodes out of sight as this could result in injury to patients / users.
- Don't plug the controller into a multi-socket outlet at the same time as other equipment. This reduces the level of safety.
- Don't allow the switched-on electrodes to come into contact with devices that conduct electric current, don't place the electrodes in the immediate vicinity of these devices.
- Don't use damaged electrodes.
- When using electrodes with a suction tube, make sure that the suction tube of electrode is connected to vacuum extractor of liquids with a possible vacuum adjustment (see p. "Safe combinations with other products"). Turn on the vacuum extractor and adjust the negative pressure from a minimum of 200 mm Hg. (26.6 kPa) to a maximum of 400 mm Hg. (53.2 kPa).
- Before carrying out operations in the depressions or cavities of patient (for example, the umbilical ring, vagina), make sure that there are no particles of flammable materials, reagents in them.
- Don't connect wet components to the controller.
- The surgical system can't handle non-conductive fluids. Don't use distilled water to irrigate electrodes. Use only sterile **conductive** solutions (e.g. saline).
- To avoid the risk of electric shock, the equipment must only be connected to a mains supply with protective earthing.
- Cautions on EMC (see section 19).

Precautions

- Before using the surgical system for the first time, carefully read all sections of Operation Manual.
- The safety and efficiency of electrosurgical procedures depends not only on the design of equipment, but also on the factors for which the user is responsible: monitoring the state of system and timely maintenance. Only persons with appropriate training and experience in electrosurgical procedures should perform procedures with this surgical system.
- It is prohibited to use electrodes to remove bone surfaces.
- It is prohibited to use the surgical radio-frequency plasma system if there is no interstitial electrically conductive medium in the area of surgical intervention or if there is no supply of saline to the electrode.
- Before performing any surgical procedures with the system, the surgeon must complete a special training course in working with electrosurgical equipment.
- Before performing electrical surgery, the physician should assess the patient's medical problems to avoid possible side effects.
- The physician must have a thorough understanding of principles and methods associated with electrosurgical procedures in order to avoid possible risks of use, as well as to prevent damage to the device and other medical products.
- If, when working with the system, it becomes necessary to use additional devices from different manufacturers, it's necessary to check their compatibility before starting the procedure.
- The electrodes are supplied sterile. Before using the electrode, check the expiration date on the packaging. Dispose of electrode if its packaging is damaged.
- The electrodes are for single use only. Don't clean, sterilize or reuse the disposable electrode.
- The electrode should only come into contact with the surgical site. Don't allow the electrode to come into contact with non-target tissue. Don't turn on the controller until the electrode is in close proximity or direct contact with the target tissue.
- When using electrodes during joint surgery, ensure that the electrode head is completely surrounded by conductive liquid (saline).
- When using electrodes with irrigation and suction, ensure continuous circulation of conductive fluid to avoid excessive heating and thermal tissue damage.
- Electrodes are not designed for cardiac radio frequency ablation.
- At the end of surgical procedure, remove the electrodes from the surgical site and place them away from metal objects.
- Don't wrap the electrode cable around metal objects; otherwise, it may result in electric shock, fire, or injury to the patient or surgical personnel.
- Don't use the electrode as a lever to increase the surgical field or gain tissue access. This can lead to bending or disconnection of electrodes and damage to the device. The electrode is not intended for mechanical tissue displacement.
- Don't allow the patient to come into contact with grounded metal objects such as the desktop frame or dashboard to avoid potential electric shock. Don't use grounding pads.
- Don't connect metal objects to an activated electrode.
- Don't use flammable anesthetics or oxidizing gases such as nitric oxide and oxygen during the procedure.
- Don't use flammable substances to clean or disinfect the generator or cables.
- High frequency electrosurgical equipment can interfere with other electronic equipment.
- The user is not allowed to disassemble the generator by himself. Repairs may only be carried out by qualified service personnel.
- Before each use, make sure that all light and sound indicators indicate normal operation of system functions. Make sure the power plug is properly connected to the generator outlet.
- When replacing the controller fuse, use the same type and class fuse that was supplied.
- The surgical system should only work with electrodes manufactured by JIANGSU BONSS MEDICAL TECHNOLOGY CO. LTD specifically for this system.

- Set and maintain the lowest power level sufficient to achieve the required procedure.
- The surgical system and electrodes should be stored away from sources of heating and direct sunlight. (see section "Storage conditions")

6. TECHNICAL CHARACTERISTICS

6.1. Electrical and technical characteristics of surgical system

Rated voltage	$\sim 230 \text{ V} \pm 10\%$
Network frequency	50 Hz
Rated current	5 A
Fuse type	RF1-20 (5A)
Power consumption of all equipment	$\leq 800 \text{ W}$
Output frequency	$100 \text{ kHz} \pm 10 \text{ kHz}; -$

Output power

Execution version	Output power in the ABLATE mode (W) $\pm 20\%$	Output power in the COAG mode (W) $\pm 20\%$	Adjusting the activation time
ARS800	380	60	-
ARS700	330	60	-
ARS600	330	60	0-900 ms, regulated


Output peak voltage

Execution version	in the ABLATE mode V_p (V)	in the COAG mode V_p (V)
ARS800	786	280
ARS700	695	280
ARS600	695	280


6.2. Electrical and technical characteristics of electrodes.



Rated voltage (all electrodes)	$\sim 500 \text{ V}$
Rated output voltage	$\sim 0 - 316 \text{ V} \pm 20\%$;
Rated output current	1A
Output frequency	$100 \text{ kHz} \pm 10 \text{ kHz};$
Type	BF

6.3. Effective output power curve when the power value is set on the display to the value "5" (half setting) and "9" (full setting) at a load of 10 Ohm - 1000 Ohm


6.4. No load power output when set to values "2" to "9"

6.5. Output power at 250 Ohm load when set to values from "2" to "9"

6.6. Classification

Type of medical device in accordance with the nomenclature classification of medical devices 236940 (Radio frequency ablation system).

Class I by type of protection against electric shock

Type BF, marked with symbol

Protection against moisture and dust:

- Shift pedals IPX8.
- IP X2 controller

Operating mode: continuous (continuous operation with intermittent load: 40 s on, 10 s off).

- - not applicable for use near flammable anesthetics and in environments with an increased oxygen content

Operation conditions

The surgical system is intended for use in the surgical departments of medical and preventive treatment institutions.

Ambient temperature	10 - 40°C
Relative humidity	30 ~ 75%
Air pressure	700 - 1060 gPa

7. INSTALLATION ORDER

7.1. Assembly diagram

The surgical radio-frequency plasma system consists of the following main parts: Controller (three versions: models ARS600, ARS700, ARS800); power cord; shift pedals AJ120 (have two modes of operation - ablation (resection) and coagulation (hemostasis of blood vessels); flow control unit and cable of flow control unit (as part of ASR 700), responsible for the flow of saline; fuses (type RF1-20 (5A), as well as various types of radio frequency plasma surgical electrodes. ARS 700 has a socket for connecting a flow control unit (used with electrodes that have an additional tube for irrigation - electrolyte delivery). The assembly diagram of surgical system is shown in (Fig. 1).



Fig. 1.
1. Controller
2. Power cord
3. Switch pedal button AJ120
4. Ablation pedal

5. Coagulation pedal
6. Shift pedal cable
7. Shift pedal cable plug
8. Cable of surgical electrode
9. Surgical electrode
10. Extraction tube of electrode (if available)
11. Irrigation tube of electrode (if available)
12. Electrode cable plug
13. Cable plug of flow control unit (for ARS 700)
14. Cable of flow control unit (for ARS 700)
15. Flow control unit (included in ARS 700)

7.2. Control panel and symbols

Front panel


1. Connector symbol for electrode plug
2. Socket for electrode plug

3. Logo of manufacturer BONSS
4. Connector for flow control unit (for ARS700)
5. Symbol of connector for flow control unit (for ARS700)
6. Indicator of connecting flow control unit (for ARS700) (red - no connection, green light - correct connection)
7. Indicator of connecting electrode (red - no connection, green light - correct connection)
8. Ablation mode control key (increase)
9. Ablation activation indicator (the indicator will light up yellow when the electrode is connected and ablation mode is activated)
10. Symbol of ablation mode
11. Power connection indicator
12. Coagulation activation indicator (the indicator will light up blue when the electrode is connected and the coagulation mode is activated).
13. Symbol of coagulation mode
14. Coagulation mode control key (increase)
15. Number of controller model
16. On/Off button (turns the AC power on and off. When turned on, the display on the controller is active. The display can remain active for 5 seconds after the power is turned off).
17. Coagulation mode control key (decrease)
18. Coagulation window on the display (shows the output coagulation level. Nominal settings will be automatically disabled when the controller is turned on with the electrode attached or when a new electrode is connected. When the controller is turned on and before the electrode is connected, the window will show "0"). There is a coagulation symbol under the window 
19. Time window on the display. (When the electrode is connected, the display window will show the working time of electrode).
20. Window of ablation on the display (shows the output level of ablation). The nominal settings will be automatically displayed when the controller is turned on with an attached electrode or when a new electrode is connected. When the controller is turned on and before connecting the electrode, the window will show "0". There is an ablation symbol under the window .
21. Ablation mode control key (decrease)
22. Light indicator for turning on. Lights up green when the controller is on.
23. Indicator for switching pedals (the indicator will light red when there are no pedals connected, green - when connected correctly).
24. Connector for the plug of switch pedal cable.
25. Connector symbol for switch pedal plug. 

Rear panel


26. Power cord connector (use only the supplied 3-wire grounded power cord).
27. Grounding symbol. 
28. Nameplate
29. Volume control for sound signals (turning the control clockwise increases the volume; counter-clockwise - decreases the volume)
30. Fuses.

7.7. Weights, dimensions, functional characteristics


Name	Appearance	Dimensions	Weight
Controller Models: ARS600 ARS700 ARS800		For all models: Length $47,0 \pm (0,4)$ cm Width $40,0 \pm (0,4)$ cm Height $13 \pm (0,2)$ cm	For all models: $9,5(-0/+1,5)$ kg
Shift pedals AJ120		Length of cable $5,0 \pm 0,1$ m Maximum width 34 cm Maximum height 6,5 cm	2200 ± 50 g
Fuse 5A		Diameter 5 mm Length 20 mm	$1 \pm 0,1$ g
Flow control		Length at the longest part $22,5 \pm 1$ cm	1350 ± 50 g

OPERATION MANUAL


OPERATION MANUAL

unit (included in ARS 700)		Width at the widest part 10.8 ± 0.5 cm Height 9.2 ± 0.5 cm				
Connecting cord of flow control unit		Length 310 ± 2 cm		160 ± 7 g		
Power cord		Length 310 ± 2 cm		270 ± 7 g		
Electrodes for ARS 600 (orthopedic surgery) Cable length of all electrodes $3 \text{ m} \pm 0.02 \text{ m}$		Irrigation	Extraction	Working head	Appliance	Weight and dimensions
1. Electrode AC301, 301Q2		NO	NO		Nucleoplasty of cervical and lumbar spine	Net weight 50 ± 5 g Electrode rod length 101 ± 5 mm Guide needle length 101 ± 3 mm Electrode rod diameter 1 ± 0.2 mm Diameter of guide needle 1.2 ± 0.2 mm Bending angle - non



OPERATION MANUAL

2. Electrode AC302, 302Q4;		NO	NO		Nucleoplasty of cervical and lumbar spine	Net weight 50 ± 5 g Electrode rod length 218 ± 5 mm Guide needle length 156 ± 5 mm Electrode rod diameter 1.1 ± 0.2 mm Diameter of guide needle 1.5 ± 0.2 mm Bending angle - none
3. Electrode MC302, 302Q5		NO	NO		Nucleoplasty of cervical and lumbar spine. Reaches the spine through the instrument channel in the endoscope	Net weight 113 ± 10 g Electrode rod length 395 ± 10 mm Electrode rod diameter 2.4 ± 0.4 mm Bending angle - none



OPERATION MANUAL

4. Electrode AC306, 306Q1;		NO	NO		Chronic tendon diseases (tennis elbow, aponeurosis, shoulder peliarthritis)	Net weight 90 ± 10 g Electrode rod length 84 ± 4 mm Electrode rod diameter 1.4 ± 0.3 mm Bending angle - $24^\circ \pm 10^\circ$
5. Electrode AC306, 306Q2		NO	NO		Arthroscopy of the hip joint Effective for focused incisions around sensitive structures in the thigh tissue. Able to eliminate sharp edges in the joint.	Net weight 90 ± 10 g Electrode rod length 130 ± 5 mm Electrode rod diameter 3.8 ± 0.5 mm Bending angle $24^\circ \pm 10^\circ$




OPERATION MANUAL

6. Electrode AC 404, 404Q1		NO	YES		Arthroscopy of the knee, shoulder, hip. Electrode for efficient ablation	Net weight 115 ± 10 g Electrode rod length 130 ± 5 mm Electrode rod diameter 3.9 ± 0.5 mm Bending angle 35° ± 10°
7. Electrode AC404, 404Q2		NO	YES		Arthroscopy of the knee, shoulder joint, meniscectomy, The shape is suitable for treating large surfaces.	Net weight 115 ± 10 g Electrode rod length 130 ± 5 mm Electrode rod diameter 3.9 ± 0.5 mm Bending angle 35° ± 10°


OPERATION MANUAL

8. Electrode AC405, 405Q2;		NO	YES		Arthroscopy of the knee, shoulder, hip joint, meniscectomy. For meniscectomy Performs resection and ablation while contouring tissue edges	Net weight 105 ± 10 g Electrode rod length 135 ± 5 g Electrode rod diameter 3.9 ± 0.5 mm Head bending angle $90^\circ \pm 10^\circ$
9. Electrode AC405, 405Q3		NO	YES		Arthroscopy of the knee, shoulder, hip joint, meniscectomy. Electrode for processing large-area joint	Net weight 105 ± 10 g Electrode rod length 135 ± 5 mm Electrode rod diameter 3.9 ± 0.5 mm Head bending angle $90^\circ \pm 10^\circ$
Electrodes for ARS 700 (ENT surgery) Cable length of all electrodes $3 \text{ m} \pm 0,02 \text{ m}$		Irrigation	Extraction	Working head	Appliance	Weight and dimensions




OPERATION MANUAL


10. Electrode AC304; A4845S		NO	NO		Uvulopalatoplasty , turbinoplasty	Net weight 109 ± 10 g Electrode rod length 117 ± 5 mm The diameter of electrode rod in the middle part is 1.6 ± 0.3 mm Bending angle - $22^\circ \pm 10^\circ$
11. Electrode AC305, A4855S		NO	NO		Uvulopalatoplasty, reduction of soft palate, uvula	Net weight 109 ± 10 g Electrode rod length 117 ± 5 mm The diameter of electrode rod in the middle part is 1.6 ± 0.3 mm Bending angle - $45^\circ \pm 10^\circ$

OPERATION MANUAL

12. Electrode AC401, A5874		YES	YES		Tonsillectomy, adenoidectomy	Net weight 160 ± 15 g Electrode rod length 120 ± 5 mm Electrode rod diameter 6.0 ± 0.5 mm Bending angle 35° ± 10°
13. Electrode AC401, A5875		YES	YES		Tonsillectomy, adenoidectomy	Net weight 160 ± 15 g Electrode rod length 110 ± 5 mm Electrode rod diameter 6.0 ± 0.5 mm Bending angle 40° ± 10°
14. Electrode AC402, A8872		YES	YES		Tonsillectomy, adenoidectomy	Net weight 160 ± 1 g Electrode rod length 120 ± 5 mm Electrode rod diameter 4.7 ± 0.3 mm Bending angle 35° ± 10°

OPERATION MANUAL


15. Electrode AC 403, A7070		YES	YES		Papillary neo- plasms of the lar- ynx, epiglottis cyst	Net weight 160 ± 15 g Electrode rod length 192 ± 8 mm Electrode rod diame- ter 3.8 ± 0.2 mm Bending angle 43° ± 10°
16. Electrode AC403, A7171;		YES	YES		Papillary neo- plasms of the lar- ynx, epiglottis cyst	Net weight 160 ± 15 g Electrode rod length 192 ± 8 mm Electrode shaft di- ameter 4.2 ± 0.2 mm Bending angle 43° ± 10°
Electrodes for ARS 800 (urology, gynecology) Cable length of all electrodes 3 m ±0,02 m		Irriga- tion	Extrac- trac- tion	Working head	Appliance	Weight and dimen- sions
17. Electrode AC307, 106Q1		NO	NO		Cutting loop elec- trode for urology and gynecology	Net weight 112 ± 10 g Electrode rod length 279 ± 5 mm Loop width in the long section 0.35 ± 0.05 mm

OPERATION MANUAL					
18. Electrode AC307, 106Q2		NO	NO		Cutting loop electrode for urology and gynecology Net weight 112 ± 10 g Electrode rod length 279 ± 5 mm Loop width in the long section 0.7 ± 0.05 mm

Plug diameter of electrodes of AC series	15±3 mm
Plug diameter of electrodes of MC series	10±3 mm

8. TYPES OF ELECTRODES

Each electrosurgical electrode has a cable (1) with plug (2) for connection to the socket of surgical system, a handle (3), an electrode rod (4), the tip of which, the head (5) is the working part. A low-temperature plasma layer is formed in the area of electrode head.


8.1. Electrodes with guide needle

This type of electrode is used with the ARS 600 controller for cold plasma nucleoplasty.


The guide needle is a thin tube that acts as a directive for the needle electrode. Under the control of intraoperative X-ray, a guide needle with an electrode is inserted into the nucleus of intervertebral disc.

8.2. Electrodes with extraction tube


The electrodes have an extraction tube and must be connected to vacuum extractor to suck out the decay products during the surgical procedure.

8.3. Electrodes with irrigation component and extraction tube

Electrodes with an additional irrigation component (irrigation system) for delivery of saline to the working head, as well as a suction tube (for suction of liquid and removed particles). (See section "Safe combinations with other products").

The flow control valve allows decreasing or increasing the delivery of saline solution to the working head. When the valve wheel is lowered down, the tube is pinched and the flow is cut off. As an electrolyte, the manufacturer recommends using a physiological solution (NaCl 0.9%) in plastic bags with a volume of 500 ml or more (for example, LS-001156 of July 10, 2019, "Sodium chloride" produced by "IST-PHARM" LLC, Russia.

9. LIST OF MATERIALS, COMING INTO CONTACT WITH AN ORGANISM

The product doesn't contain medicines
The product doesn't contain materials of animal or human origin.

No.	Name of parts and assemblies, execution options		Type and grade of material.
1	Controller	Metal case	"XURI" galvanized iron, manufactured by Suzhou XuriGuantong Metal Products Co Ltd
		Control panel	ABS plastic «HUANXIN», manufactured by Suzhou HuanXin Precision Mold Co Ltd
2	Flow control unit	Metal case	"XURI" galvanized iron, manufactured by Suzhou XuriGuantong Metal Products Co Ltd
		Cable	«TAOSHI» silicon rubber, manufactured by Dongguan RiZhenShangQin Electrocn Materials Co Ltd
		Plug	Polytetrafluoroethylene «DAJIN»,

			manufactured by Dongguan RiZhenShangQin Electrochemical Materials Co Ltd
3	Electrodes	Tip of electrode (head)	Ceramics «ZHAOHE», manufactured by Shenzhen FangTai New Material Technology Co Ltd Stainless steel «CHENGZHIDE 304», manufactured by Shenzhen Chengzhide Precision Tube Co Ltd, PVDF (polyvinylidene fluoride) «CHANGYUAN» brand, manufactured by Shanghai Changyuan Electronic Materials Co Ltd,
		Rod of electrode (for uncoated electrodes)	Stainless steel «CHENGZHIDE 304», manufactured by Shenzhen Chengzhide Precision Tube Co Ltd
		Coating of rod (for coated electrodes)	PVDF (polyvinylidene fluoride, polytetrafluoroethylene), CHANGYUAN brand, manufactured by Shanghai Changyuan Electronic Materials Co Ltd Coloring agent «BASF black» brand, manufactured by Yangzi Shihua BASF Limited Liability Company
		Handle (holder) and Plastic plug for electrode cord	ABS plastic, «HUANXIN», manufactured by Suzhou HuanXin Precision Mold Co Ltd Coloring agent: CHIMEI blue CHIMEI black, Manufacturer - Zhenjiang ChiMei Chemical Co Ltd
		Guide needle, Tube of guide needle	Stainless steel YingFeiNiAo, manufactured by Jiangzhou Medical Device Company
		Handle of guide needle	- Polycarbonate (PC), JinYi brand, Jiangzhou Medical Device Company
		Yellow stopping device	Thermoplastic plastics TPU, JinYi brand, coloring agent «BASF yellow» brand,

		Green stopping device	manufactured by Jiangzhou Medical Device Company Thermoplastic plastics TPU, JinYi brand, coloring agent «BASF green» brand, manufactured by Jiangzhou Medical Device Company
		Трубка для эвакуации (suction)	PVC (polyvinyl chloride) «YINUO» brand, manufactured by Nantong YiNuo Precision Plastic Tube Co Ltd
		Cable of electrode,	Silicon rubber, «TAOSHI» brand, manufactured by Dongguan RiZhenShangQin Electrocn Materials Co Ltd
		Plug of electrode made of PTFE	PVDF (polyvinylidene fluoride), «DAJIN» brand, manufactured by Dongguan RiZhenShangQin Electrocn Materials Co Ltd
4	Power cord		Silicon rubber, «TAOSHI» brand, manufactured by Dongguan RiZhenShangQin Electrocn Materials Co Ltd Coloring agent CHIMEI black, Manufacturer - Zhenjiang ChiMei Chemical Co Ltd

10. SAFE COMBINATIONS WITH OTHER PRODUCTS

1. When using electrodes with an extraction tube:

1.1 (electrodes included in the execution version ARS 600:

- Electrode AC 404, 404Q1,
- Electrode AC404, 404Q2,
- Electrode AC405, 405Q2,
- Electrode AC405, 405Q3,

it's necessary to use them together by vacuum liquid extractor with possible adjustment of negative pressure from a minimum of 200 mm Hg. (26.6 kPa) to a maximum of 400 mm Hg. (53.2 kPa).

For instance,

- a pump for irrigation (as part of medical device "Video endoscopic system for carrying out arthroscopic operations with accessories" manufactured by "ConMed Corporation", USA, RU No. FSZ 2010/07098 dated October 19, 2019.
- version of "Vacus 7032" medical device "Vacuum extractor "Vacus" according to TU 9444-005-74487176-2010" manufactured by "DIXION" LLC, Russia, RU No. FSR 2010/08658 dated October 07, 2019).

Attention! The flow control unit included in the ARS 700 execution version is not compatible with the ARS 600.

1.2. electrodes included in the execution version ARS 700:

- Electrode AC401, A5874,

- Electrode AC401, A5875,
- Electrode AC402, A8872,
- Electrode AC403, A7070,
- Electrode AC403, A7171

it's necessary to use them together with medical vacuum extractor of liquids with possible vacuum adjustment from 0.02 MPa - 0.09 MPa.

For example, the execution version "Vacus 7032", the medical product "Vacuum extractor "Vacus" according to TU 9444-005-74487176-2010" manufactured by "DIXION" LLC, Russia, RU No. FSR 2010/08658 dated October 07, 2019).

2. When using an electrode requiring insertion through the instrument channel of endoscope:

2.1. Electrode as part of execution version ARS 600

- Electrode MC302, 302Q5

it's necessary to use it together with endoscopic medical equipment, which includes an endovideo camera, a light source, an endoscope (for example, a medical device "Video endoscopic system for carrying out arthroscopic operations with accessories" manufactured by "ConMed Corporation", USA, RU No. FSZ 2010 / 07098 dated October 19, 2019).

3. When using all other electrodes that are part of surgical system options, combinations with other products are not required.

11. SOFTWARE

The program is part of medical device. The consumer receives a medical device with an installed program.

11.1. System configuration

Microcontrollers	Auxiliary Atmel atmeg8	Main Atmel atmeg32
Maximum frequency	8 MHz	8 MHz
Random access memory (RAM)	1 Kbyte	2 Kbyte
Flash memory	8 Kbyte	32 Kbyte
Nucleus	AVR	AVR

The software used in the medical device is installed by the manufacturer on the flash memory of Atmel atmeg8 (auxiliary) and Atmel atmeg32 (main) microcontrollers

11.2. Firmware elements

Identifier of program	ARSdriverV01M for the main program ARSdriverV01S for the auxiliary program
Version of operation system	Not used
Version of program	Not less than ARSdriverV01.01M Not less than ARSdriverV01.01S
Hardware version of program	BONSSMAIN V04
Date	June 15, 2017
Safety class	B
Total size of program	Less than 30 Kbyte (ARSdriverV01.01M – less than 20 Kbyte, ARSdriverV01.01S – less than 10 Kbyte)
Search time of program	Not applicable
Response time of program	5 s

11.3. Information on errors and statuses

Malfunctions without error code reporting

Error and status information is shown on the information display. It is reflected as red and green LEDs, and there are operational audible signals and alarms.

Description of error	Description of status data/solution of problem
When turning on the surgical system, no information is displayed on the control panel	No power. 1. Check the supply voltage AC 230 V $\pm$ 10%; 2. Check if the power cable is damaged; 3. Check if the fuse is damaged.
After the system is activated, no plasma layer forms on the tip of surgical electrode	1. Check the electrode irrigation system. 2. Confirm the successful connection of electrode, cables and foot switching pedal. 3. Perform a test using a new cable, new electrode, or new foot pedal to ensure they are working properly. If the components have been successfully tested, contact an Authorized Service Center.
Digital display «CE»	Indicates no connection with the electrode. Check the correct connection of electrode plug. Try to connect again.
Digital display «00»	Waiting state for electrode connection. Connect the electrode to controller connector
The electrode indicator doesn't turn green when the electrode is connected	Electrode recognition error. Check if the connection is correct. In case of repeated error, replace the electrode.
During operation, the electrode indicator lights up red. The digital display shows "00". Intermittent alarm sounds	Controller is defective / electrode is defective. Stop work. Turn off the controller. Replace electrode. If the signals are repeated, contact the service center.

Replacement of fuse:

- 1) Turn off the power and disconnect it from the main power supply, insert the slotted screw into the bolt on the fuse cover, rotate 90 degrees counterclockwise, and then the fuse cover can be removed.
- 2) Remove the fuse cover and the damaged fuse.
- 3) Insert a new fuse with the same parameters into the fuse body, making sure it is securely locked.
- 4) Insert the cover into the cover slot, insert the slotted screw into the bolt on the fuse cover, and rotate 90 degrees clockwise. Make sure the fuse cover is securely locked.

12. WORKING AND ALARM SIGNALS**Characteristics of sound and light signals**

12.1. Working signals: When the system is turned on, a single continuous signal is heard; the operation indicator light is on and the system starts generating radio frequency energy.

Normal ablation signal with pedal depressed - constant signal (600 $\pm$ 50) Hz

Normal coagulation signal with pedal depressed - constant (1200 $\pm$ 50) Hz

12.2. Alarm signal:

intermittent double signal (intermittent double signal at normal ablation signal frequency for 0.25 seconds, then coagulation signal for 0.25 seconds until start event disappears), operation indicator light is off, no radio frequency energy output.

13. CLEANING, DISINFECTION AND MAINTENANCE

All components of surgical system (with the exception of electrosurgical electrodes) are supplied non-sterile. Don't sterilize during commissioning and subsequent use.

Maintenance of surgical system consists of weekly cleaning and disinfection of controller surfaces, AJ120 switch pedals, and flow control unit.

Before first use and weekly, the controller, AJ120 pedals, and flow control unit should be cleaned with a soft cloth dampened with 0.05% chlorhexidine digluconate solution. Don't spill disinfectant on the surface of controller.

Other components of product don't need to be disinfected.

Electrodes

All surgical electrodes are for single use only and are supplied sterile in sterile packaging. Sterilization method is ethylene oxide. (see "EO Sterilization Validation Report" prepared in accordance with ISO 11135: 2014, ensuring a sterility assurance level (SAL) of at least 10^{-6}).

Reuse and re-sterilization of electrodes is strictly prohibited. The manufacturer is not liable for accidents resulting from the reuse and re-sterilization of surgical electrodes.

If soft tissue particles adhere to the electrode tip during use, use a soft cloth soaked in 96% alcohol solution to clean the electrode tip. The electrode is cleaned when the controller is off. Don't use sharp tools to clean the electrode tip. Don't use or continue to use the electrode if the insulation of electrode rod is damaged.

Repair of controller

Attention! Only qualified technicians, trained and authorized by the manufacturer, can arrange for controller repair. Failure to comply with this requirement will void the client's right to receive any additional service from the manufacturer.

14. WARRANTY OBLIGATIONS

1) The manufacturer is responsible for the quality and compliance of all components of surgical technical documentation system.

2) The manufacturer is not liable under the following circumstances:

A. The presence of defects or problems arising from the violation of conditions of transportation upon delivery of product to the original purchaser.

B. Injury to patients / users and property damage caused by improper installation, use, repair or storage of product, as well as the user's failure to comply with the warnings and cautions provided by the manufacturer.

D. The customer or user doesn't comply with a written request from the manufacturer related to performance, improvement or modification of product.

E. Injury to patients / users and property damage caused by product changes or modifications that are not controlled by the manufacturer.

F. Patient / user injury and property damage caused by customer or user negligence, recklessness, or willful violation of product requirements.

15. CONDITIONS OF STORAGE AND TRANSPORTATION

Ambient temperature	-40-70 °C
Relative humidity	10~100 %
Atmospheric pressure	500-1060 gPa

Store this device in a well-ventilated area. The air must not contain aggressive gases. Storage height above floor level and distance from wall must be more than 30 cm.

Maximum stacking during storage and transportation in transport packaging in height - 6 boxes.

In case of long-term storage, the surgical system is recommended to be put into operation for a short time to check its serviceability at least 1 time per year.

Pack the surgical system before shipping to protect it from strong shocks, crushing, or dropping. The surgical system can be transported by car, train-house, ship and air transport.



16. MARKING AND PACKING

16.1. Transport packaging of surgical system

The surgical system is packed in boxes made of five-layer corrugated cardboard with dimensions of 53 x 50 x 37 (±0,5) cm. Internal foam and cardboard inserts provide the system and its components with safe transportation.

A packing list is attached to the box with information about the name of product, the name and address of manufacturer, information about the buyer and his address.

The following symbols are printed on the outside of transport packaging:


	This side up.	Indicates the correct vertical position of the cargo
	Protect from moisture.	The need to protect cargo from moisture exposure
	Fragile. Take care.	Careful handling of cargo
	Limit on the number of tiers in stack.	The maximum number of identical cargos that can be stacked one on top of the other, where 6 is the limit
	Avoid exposure to sunlight	Indicates that the medical device must be protected from exposure to sunlight

16.2. Transport packaging of electrodes


Individually packed electrodes in an amount of at least 10 pieces are placed in corrugated cardboard boxes. Box size 43.5x33.5x14.5 (± 0.5) cm. The maximum weight of transport box with 10 electrodes is 5.5 kg.

The transport packaging for electrodes carries the same symbols and information as the packaging for the surgical system.

16.3. Marking of controller

There is a label on the back of controller panel containing the following information:

The “Use only grounded receptacles” symbol is located on the rear panel of controller.

- Product name and execution version of controller (ARS 700)
- Rated supply voltage
- Rated current of the controller fuse
- Voltage frequency
- Rated output power at 250 Ohm load)
- Output frequency
- Manufacturer's name and production address
- Authorized representative in the EU
- Serial number
- Symbol "Use before"
- Date of manufacture
- Symbol “Non-ionizing radiation!”
- Symbol "Caution! Please refer to the instructions for use”.
- Symbol "Don't dispose of as household waste"
- Symbol "Refer to instructions for use"
- CE mark and registration number. The product complies with the essential requirements of Medical Device Directive 93/42 / EEC
- Degree of protection against moisture and dust IP

16.4. Label layouts for RF controllers.

On the rear side of controller panel, a label in Russian must be applied, containing the following information:

Execution version ARS 600

Surgical radio frequency plasma system ARS

ARS 600	RU No. ____	
~ 230 V	5 A	50 Hz
Output	330 W /250	100 kHz
	Ohm	
Operating mode	40 s on,	10 s off.

JIANGSU BONSS

Medical Technology Co.Ltd.
Building #7, No.898, China Medical City Avenue, Hailing District,
Taizhou City, 225316, Jiangsu P.R. China

SN

IP X2

Class I

YEAR.MONTH

YEAR.MONTH

CE 0197

Execution version ARS 700

Surgical radio frequency plasma system ARS

ARS 700	RU No. ____	
~ 230 V	5 A	50 Hz
Output	330 W /250	100 kHz
	Ohm	
Operating mode	40 s on,	10 s off.

JIANGSU BONSS

Medical Technology Co.Ltd.
Building #7, No.898, China Medical City Avenue, Hailing District,
Taizhou City, 225316, Jiangsu P.R. China

SN

IP X2

Class I

YEAR.MONTH

YEAR.MONTH

CE 0197

Execution version ARS 800

Surgical radio frequency plasma system ARS

ARS 800	PY No ____	
~ 230 V	5 A	50 Hz
Output	380 W /250	100 kHz
	Ohm	
Operating mode	40 s on,	10 s off.

JIANGSU BONSS

Medical Technology Co.Ltd.
Building #7, No.898, China Medical City Avenue, Hailing District,
Taizhou City, 225316, Jiangsu P.R. China

SN

IP X2

Class I

YEAR.MONTH

YEAR.MONTH

CE 0197

The symbol of using protective earthed socket  is located on the rear panel of controller.

16. 5. Consumer individual packaging of electrodes

The consumer packaging of electrodes has four types of sizes

37.1x10.1x2.75 ± 0,3 cm

60.6x16.1x3.7 ± 0,3 cm

42.6x13.7x3.1 ± 0,3 cm

46.1x16.1x3.3 ± 0,3 cm

Each consumer package contains a short instruction. Weight of instructions is 12 g.

The consumer packaging of electrodes contains the following information:

- Name «Radio frequency plasma surgical electrode»
- Rated supply voltage
- Rated current of the fuse,
- Voltage frequency
- Rated output voltage
- Rated output current
- Output frequency
- Manufacturer's name and production address
- Authorized representative in the EU
- CE mark and identification number of registration body. The product complies with the essential requirements of Medical Device Directive 93/42 / EEC
- Symbol "Sterilized with ethylene oxide"
- Symbol "Don't use with damaged packaging"
- Symbol "Prohibition of reuse"
- Symbol "Non-ionizing radiation!"
- Symbol "Caution! Please refer to the instructions for use."
- Symbol "Don't dispose of as household waste"
- Symbol "Refer to instructions for use"
- Model of electrode
- Catalog number of electrode
- Catalog name of electrode
- Date of manufacture (sterilization)
- Expiry date (Use before)
- Batch number (LOT)
- Made in China

16. 6. Consumer packaging of electrodes

Name of electrode	Gross weight (individually packed) (±10% g)	Weight in a ster- ile blister (±10%)	Dimensions of individual packaging LxWxH (±0,3 cm)
1. Electrode AC301, 301Q2	163	99	37.1x10.1x2.75
2. Electrode AC302, 302Q4;	167	100	37.1x10.1x2.75
3. Electrode MC302, 302Q5	364	243	60.6x16.1x3.7

Electrode AC306, 306Q2	250	172	42.6x13.7x3.1
6. Electrode AC 404, 404Q1	272	188	42.6x13.7x3.1
7. Electrode AC404, 404Q2	255	180	42.6x13.7x3.1
8. Electrode AC405, 405Q2;	261	179	42.6x13.7x3.1
9. Electrode AC405, 405Q3	255	179	42.6x13.7x3.1
Electrodes for ARS 700			
10. Electrode AC304; A4845S	266	182	42.6x13.7x3.1
11. Electrode AC305, A4855S	266	182	42.6x13.7x3.1
12. Electrode AC401, A5874	315	238	42.6x13.7x3.1
13. Electrode AC401, A5875	317	235	42.6x13.7x3.1
14. Electrode AC402, A8872	310	234	42.6x13.7x3.1
15. Electrode AC 403, A7070	380	284	46.1x16.1x3.3
16. Electrode AC403, A7171;	352	260	46.1x16.1x3.3
Electrodes for ARS 800			
17. Electrode AC307, 106Q1	152	93	37.1x10.1x2.75
18. Electrode AC307, 106Q2	155	95	37.1x10.1x2.75

21.7. Sterile packaging of electrodes

The packaging of electrode is a transparent plastic PETG blister for medical use and Tyvek 1059B cover.

Residual amount of standard sample (SS): not more than 10 µg / g complies with the control standard: ISO 10993-7-2008 "Assessment of biological effect of medical devices - Part 7: Residual ethylene oxide content after sterilization".

Characteristics of sterile packaging	Cover of sterile blister TYVEK 1059B 1*	Cover of sterile blister TYVEK 1059B 2*	Cover of sterile blister TYVEK 1059B 3*	Cover of sterile blister TYVEK 1059B 4*	Sterile blister made of plastics PETG 1*	Sterile blister made of plastics PETG 2*	Sterile blister made of plastics PETG 3*	Sterile blister made of plastics PETG 4*
Dimensions (L x W)	37x10 cm (+0/-1)%	42.5x13.5 cm (+0/-1)%	60.5x16 cm (+0/-1)%	46x16 cm (+0/-1)%	37x10 cm (+0/-1)%	42.5x13.5 cm (+0/-1)%	60.5x16 cm (+0/-1)%	46x16 cm (+0/-1)%
Width of thermal junction seam	0,9 cm ± 0,1 cm				0,9 cm ± 0,1 cm			
Density	64,4 g/m ³ (±5%)				1,27 g/cm ³ (±5%)			
Air permeability	671 ml/min (±5%)				-			
Tear resistance	1055 kPa (±5%)				54 MPa (±5%)			
Bursting strength	72,8 H/cm ² (±5%)				360 H/cm ² (±5%)			
Average thickness value	163 μm				1,4 mm			
Thickness deviation	±22 μm				±0,05 mm			
Tensile strength	14,5 MPa (±5%)				-			
Bonding strength	1,5N/15 mm (±5%)							
Note	1* - for 1. Electrode AC301, 301Q2 2. Electrode AC302, 302Q4 17. Electrode AC307, 106Q1 18. Electrode AC307, 106Q2 2* for							

erile
ster
e of
cs
7

3. Electrode MC302, 302Q5

3* for

4. Electrode AC306, 306Q1

5. Electrode AC306, 306Q2

6. Electrode AC 404, 404Q1

7. Electrode AC404, 404Q2

8. Electrode AC405, 405Q2

9. Electrode AC405, 405Q3

10. Electrode AC304, A4845S

11. Electrode AC305, A4855S

12. Electrode AC401, A5874

13. Electrode AC401, A5875

14. Electrode AC402, A8872

4* for

15. Electrode AC 403, A7070


16. Electrode AC403, A7171

21.8. Label layout of electrodes for RF

On the Tyvek lid, a label is applied in Russian, containing information:

Surgical radio frequency plasma electrode
AC404, 404Q2
~230 V 5A 50 Hz
AC ~ 0-316 V 1A 100 kHz

 **JIANGSU BONSS**
Medical Technology Co.Ltd.
Building #7, No.898, China Medical City Avenue,
Hailing District, Taizhou City, 225316, Jiangsu P.R.
China

 0197   


 YEAR, MONTH, DAY  YEAR, MONTH, DAY

 XXXXXX **Made in China**

- Name and type of electrode
- Rated supply voltage
- Rated current of controller fuse
- Voltage frequency
- Rated output voltage
- Rated output current
- Output frequency
- Manufacturer's name and production address
- CE mark and identification number of registration body. The product meets the essential requirements of Medical Device Directive 93/42 / EEC
- Symbol "Sterilized with ethylene oxide"
- Symbol "Don't use with damaged packaging"
- Symbol "Non-toxic"
- Symbol "Prohibition of reuse"
- Symbol "Non-ionizing radiation!"
- Symbol "Caution! Please refer to the instructions for use."
- Symbol "Don't dispose of as household waste"
- Symbol "Refer to instructions for use"
- Symbol "Non-pyrogenic"
- Date of manufacture (sterilization)
- Expiry date (Use before)
- Batch number (LOT)
- Made in China

Symbols and designations are shown in the following table:

	Manufacturer	
	Authorized representative in the EU	
	Use before	
	Date of manufacture (sterilization)	
	Batch number	
	Serial number	
	Sterilized with ethylene oxide	
	Prohibition of reuse	
	Pyrogenically	
	Don't use with damaged packaging	
	Non-toxic	
	Product type BF	

		
	Non-ionizing ionization	
	Don't dispose of as household waste	
	Refer to the operating documentation	
	Attention! Refer to instructions for use"	
CE 0197	CE mark and identification number of registration body	

17. OPERATION

17.1. Unpacking, installing and checking

1. Unpacking:

Unpack the box. Make sure that the surgical system has been delivered complete and in good working order. If you find any damage, contact an authorized representative in the Russian Federation to resolve the dispute.

17.2. Assembly and functioning check

2.1 Before the first use, clean the surfaces of controller, the switch pedal, and flow control unit with a soft cloth dampened with a 0.05% chlorhexidine digluconate solution.

2.2. Plug one end of the power cord into the back of generator, and plug the other end into an electrical outlet.

2.3. Press the power switch on the front of controller to toggle to the | (ON) position and wait 5 seconds.

2.4. Connect the control pedals to the front of controller, the pedal indicator light turns green.

2.6. Remove the sterile electrode from the package. Connect the electrode to the connector on the front of controller. The electrode indicator light turns green.

2.7. When the electrode is connected, the controller recognizes the electrode and sets automatically the value of ablation and coagulation power.

2.8. Prepare a glass of normal saline solution.

2.9. Use the automatically set power value shown on the display.

2.10. Press the yellow ablation pedal and then place the electrode tip in saline. If the system is working properly, orange plasma radiation is generated around the electrode tip, accompanied by a hissing sound.

2.11. Press the blue coagulation control pedal and then place the electrode tip in saline. If the system is working properly, bubbles will form around the electrode tip, accompanied by a hissing sound.

If the system is not working properly, contact the service department of authorized representative in the Russian Federation.

17.3. Instructions during operation:

The system is assembled according to fig. 1.

1. Select the appropriate electrode type. Connect it to the controller.

1.2. For electrodes with an irrigation system:

- for ARS 700 execution version, connect the electrode irrigation tube to the flow control unit, (see p. Installation of flow control unit).

- for ARS600 and ARS 800 execution versions it's necessary to connect the irrigation tube to irrigation pump (see the section "Safe combinations with other products).

1.3. For electrodes with an extraction tube, it's necessary to connect to vacuum extractor (see the section "Safe combinations with other products").

2. During the operation, the surgeon can adjust the appropriate power parameters for the ablation and coagulation modes.

3. Position the working electrode head near the surgical site, turn on compatible devices, and press the switch pedal to activate the device.

Check the working electrode head periodically to make sure the tip is not damaged and the electrode is working as instructed. If the electrode tip is intact, discontinue use and check the surgical site accordingly. Overuse and / or use above the default preset value can lead to electrode breakage and shorten the life of device.

2. If any residual material adhered to the electrode head during single patient surgery, wipe the head with gauze and rubbing alcohol. Any other items are not recommended to be used to remove residual materials.

17.4. Instructions after operation:

1. System shutdown: press the power switch and wait 5 seconds until all generator indicators turn off.

2. Remove the plug from the socket.

3. Disconnect the electrode cable from the controller.

4. Dispose of electrodes.

17.5. Installation of flow control unit (for ARS 700)

1. Front side of the flow controller.

2. Fasten the control box to a rigid surface with the screw on the back of flow controller.

3. Connect one plug of cable to the connector on the front of controller and the other end of cable to the connector on the rear of flow control unit. Make sure the cable is connected in the correct direction. After successful connection, turn on the power switch of flow control unit.

4. When turned on, the power indicator on the front of flow control unit is solid green and the flow control valve opens to insert the saline irrigation tubing.

5. Insert the saline tubing into the flow control valve and turn off the switch on the flow control unit.

6. Open the valve on the flow control tubing and saline starts to flow. During the procedure, the flow controller can be activated at the same time as the controller using the switch pedals.

18. ENVIRONMENTAL PROTECTION

The surgical system (controllers, power cord, cable of flow control unit is subject to the European RoHS directive (prohibition of the use of certain hazardous substances in the Electrical and Electronic Equipment Directive).
Related information

Name of detail	Toxic and harmful substances or elements					
	Pb	Hg	Cd	CrVI	PBB	PBDE
Mother board	*					
Power cord	*					
Cable of flow control unit	*					

- components in all toxic or harmful substances with the same material quality below the specified limit.
- * components of at least one type of toxic or harmful substances with a similar material quality above the specified concentration limit.

This device complies with EU Directive 2011/65/EC (provisions on restrictions on the use of certain hazardous substances).

This equipment must be returned by the user to an authorized company for disposal or recycling.

19. ELECTROMAGNETIC COMPATIBILITY (EMC)

This equipment has been tested and found to comply with the limits for:
EN 60601-1-2: 2007 Medical electrical equipment. Part 1-2: General requirements for safety and essential performance assurance. Standard: Electromagnetic compatibility. Requirements and tests.

- 1) These limits are designed to provide reasonable protection against harmful electromagnetic or other interference in most installations. However, there is no guarantee that interference will not occur in a particular installation. If this equipment causes harmful electromagnetic interference or other interference, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
- A. Change direction or relocate the ARS controller
 - B. Increase the distance between the ARS controller and related equipment.
 - C. Connect non-surgical equipment to an outlet on a circuit different from that to which the system is connected.
 - D. Consult the dealer or authorized technical personnel for help.

TABLES AND DECLARATIONS REGARDING ELECTROMAGNETIC COMPATIBILITY.

Table 2: Guidance and declaration of manufacturer - electromagnetic emissions
Surgical radio frequency plasma system ARS.
Guidance and declaration of manufacturer - electromagnetic compatibility
Classification

Regulatory document	Group	Class
GOST R 51318.11	1	A
GOST 30804.3.2		A
GOST 30804.4.11		2
GOST R IEC 60601-1-2	ME PRODUCT or ME SYSTEM, CONNECTED with PATIENT	

The system is intended for use in the electromagnetic environment specified below.

Guidance and declaration of manufacturer - electromagnetic compatibility			
Stability test	EN 60601 Test level	Compliance level	Compliance level Electromagnetic environment
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors: ceramic tiles Humidity: 55 %
Fast electrical transient processes or bursts IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input / output lines	± 6 kV contact ± 8 kV air	The quality of electrical network should be that of a typical commercial or hospital environment.
Line voltage surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	The quality of electrical network should be that of a typical commercial or hospital environment.
Voltage losses, short interruptions and voltage fluctuations on input power lines IEC 61000-4-11	<5% UT (> 95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (> 95% dip in UT) for 5 sec	<5% UT (> 95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (> 95% dip in UT) for 5 sec	The quality of electrical network should be that of a typical commercial or hospital environment. If the user requires continued operation during a power outage, it is recommended that the amplifier of surgical system be connected to an uninterruptible power supply or battery.
Magnetic field frequency (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	If malfunctions occur, it may be necessary to place the Surgical System further from sources of magnetic power frequency fields or install magnetic shielding. The power frequency magnetic field should be measured at the intended installation location to ensure that it is sufficient.

Guidance and declaration of manufacturer – electromagnetic compatibility

Stability test	EN 60601 Test level	Compliance level	Compliance level Electromagnetic environment
Conducted RF IEC 61000-4-6	3 Vrms from 150 kHz to 80 MHz out of ISM range	3 Vrms from 150 kHz to 80 MHz	Portable and mobile radio frequency equipment should be used away from the surgical radio frequency plasma ARS system, including cables, than the recommended distance calculated from the equation applicable to the frequency of transmitter. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V / m from 80 MHz to 2,5 HHz	3 V / m from 80 MHz to 2,5 HHz	

Recommended distances between portable and mobile radio frequency equipment and the surgical radio frequency plasma ARS system.

Guidance and declaration of manufacturer – electromagnetic compatibility			
The surgical system is intended for use in electromagnetic environment in which radiated radio frequency disturbances are controlled. The user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile radio frequency communications equipment (transmitters) and the surgical radio frequency plasma ARS system, as recommended below, in accordance with the maximum output power of communications equipment.			
Rated maximum output power of transmitter	Separation distance by transmitter frequency (m)		
	from 150 kHz to 80 MHz $d=1.2\sqrt{P}$	from 80 MHz to 800 MHz $d=1.2\sqrt{P}$	from 800 MHz to 2,5 HHz $d=1.2\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters with a maximum output power not listed above, the recommended distance d in meters (m) can be estimated using the equation applicable to the frequency of transmitter, where P is the maximum output power of transmitter in watts (W). NOTE 1. At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2. These recommendations may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

20. SERVICE LIFE, WARRANTY OBLIGATIONS

The shelf life of radio frequency plasma system is 6 years from the date of manufacture.

The shelf life of sterile electrodes is 3 years from the date of sterilization.

The warranty storage / operation period for the surgical system is 12 months from the date of delivery to the buyer, the warranty for the switch pedals is 180 days.

Problems arising with the components of surgical system during the warranty period, associated with defective materials or equipment malfunction, the manufacturer guarantees to solve by replacing the defective equipment.

Obligations of manufacturer

Any components that don't meet the requirements of technical documentation of the surgical radio frequency plasma system, the manufacturer undertakes to repair, replace or return free of charge.

The manufacturer assumes responsibility for loss of life and property damage caused by product defects, except for the following situations:

- A. Damage resulting from violation of conditions of transportation and storage of surgical system.
- B. Non-observance by the user of precautions, warnings indicated in the operation manual.
- C. Injuries to patients / users and property damage caused by non-compliance with the conditions of transportation, storage, installation, operation, as well as repairs outside an authorized service center.
- D. Patient / user injury and property damage caused by user negligence, recklessness or willful violation of product requirements.

21. RECLAMATIONS

Complaints and claims about the quality of surgical radio frequency plasma system should be sent to the authorized representative in the Russian Federation. "SMK-Trade" Limited Liability Company ("SMK-Trade" LLC),

Russia, 214010, Smolensk region, city of Smolensk, Sedov st., 42

Tel. +7 905 163 58 60

E-mail: aleksey.zhdanyuk@yandex.ru

22. DISPOSAL

The product can't be disposed of as household waste. All components of system must be sent to the appropriate services for disposal and recycling. The device contains an integrated circuit and must be disposed of in accordance with local regulations and policies for the rejection of electronic devices.

Disposable electrodes after use are subject to disposal as medical epidemiological hazardous waste of class B (in accordance with the requirements of SanPiN 2.1.7.2790-10).

23. INTERNATIONAL STANDARDS

MDD 93/42/EEC:1993	Medical devices directive
2007/47/EC:2007	Medical devices directive
EN ISO 13485:2012	Quality system - medical instruments
EN ISO 14971:2012	Medical devices - Applying risk management to medical devices
EN 1041:2008+A1:2013	Information supplied by manufacturer for medical devices
EN ISO 10993-1:2009/AC:2010	Biological evaluation of medical devices - Part 1: Evaluation and Testing
EN ISO 10993-5:2009	Biological evaluation of medical devices - Part 5: Cytotoxicity tests in vitro

EN ISO 10993-10:2010	Biological evaluation of medical devices – Part 10: Testing on skin irritation and allergy
EN 60601-1:2006/AC:2010	Medical electrical equipment - Part 1: General safety requirements
EN 60601-1-2:2007	Medical electrical equipment - Part 1-2: General safety requirements - Parallel standard: Electromagnetic compatibility - Requirements and tests
EN 60601-2-2:2009+A11:2011	Medical electrical equipment - Part 2-2 Particular safety requirements for high-frequency electrosurgical devices
EN ISO 15223-1:2012	Medical devices - Symbols used for marking on medical devices, labels and accompanying documents - Part 1. Basic requirements
BS EN 868-1:1997	Materials and systems for packaging medical devices intended for sterilization. General requirements and test methods
BS EN 868-2:2009	Packaging of medical devices subject to terminal sterilization. Packaging for sterilization. Requirements and tests
EN ISO 14155:2011	Clinical trials of medical devices for human use - Good Clinical Practice
EN 62366:2008	Medical devices – Design of medical products taking into account serviceability
BS EN ISO 14644-1:2010	Cleanrooms and associated controlled environments - Part 1: Classification of air cleanliness
EN ISO 14698-1:2003	Cleanrooms and associated controlled environments: Biocontamination control. General principles and methods.
EN ISO 11135-1:2007	Sterilization of medical devices: Ethylene oxide - Part 1: Requirements for the development, validation and routine control of sterilization process for medical devices
EN ISO 11607-1:2009	Packaging of terminally sterilized medical devices - Part 1: Requirements for materials, sterilization barrier systems and packaging systems
EN ISO 11607-2:2006	Packaging of terminally sterilized medical devices - Part 2: Requirements for validation of forming, sealing and assembly processes
EN ISO 11737-1: 2006	Sterilization of medical devices - Microbiological techniques - Part 1: Evaluation of population of microorganisms on products
EN 60601-1-6:2010	Medical electrical equipment. Part 1-6: General safety requirements with regard to essential functional characteristics. Parallel standard: Serviceability
EN 62304: 2006/AC:2008	Software of medical devices - Life cycle processes of software means

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