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102
Jut 59673

Registered No. 2021 - 731

NOTARIAL CERTIFICATE

CHUNG SOL LAW AND NOTARY OFFICE

365 Sinwol-ro, Yangcheon-gu, Seoul, Korea
(Sinjeong-dong, Chongsol Building 2F)

Statement

ProteomeTech Inc. declares that we issued Operational manual as attached.

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
ProteomeTech Inc.

Date: June 8th, 2021


LIMITED LIABILITY COMPANY
"Bio Focus"

AUTOMATED DEVICE "Q-PROCESSOR" FOR "PROTIA" PANELS PROCESSING FOR IN VITRO DIAGNOSTICS

Operational manual

26.60.12.002.15580960 RE

Moscow
2021

**DATA FOR PRODUCT IDENTIFICATION**

Type and designation: *Automated device "Q-PROCESSOR" for "PROTIA" panels processing for in vitro diagnostics*

Model: _____

Serial number: 6519031006

Manufacturing date: 03.2019

Name and address

- manufacturer Company "ProteomeTech Inc."
A-702, 401 Yangcheon-Ro, Gangseo-Gu, Seoul 07528, Republic of Korea
+ 82-2-6968-5277 / +82-2-6968-5279

- seller (supplier) LLC "Bio Focus"
109387, Russia, Moscow, Tikhaya str., 28, office 216

Phone fax: +7 (916) 900-65-84

Website, e-mail:

general issues: www.biofocus.ru, www.proteometech.com

technical issues: protia@proteometech.com, info@biofocus.ru

service center: biofocus@yandex.ru

Device complies with the current regulatory and technical documentation, and is recognized as fit for use

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CONFORMITY INFORMATION

State registration certificate:	110-81-45244
Issued by:	30.03.2000
Device meets the requirements of regulatory documents:	GOST R 15.013-2016, GOST ISO 13485-2017, GOST 20790-93, GOST R 50444-2020, GOST IEC 61010-2-101-2013, GOST IEC 60065-2013, GOST IEC 60950-1-2014, GOST R IEC 60601-6-2014, GOST R 55991.7-2014, GOST R IEC 62366-2013, GOST 30324.0-95 / GOST R 50267.0-92 and GOST 12.2.091-2012

ATTENTION OF CONSUMERS!

Thank you behind acquisition Devices, manufactured our now!

We ask you to carefully study this Operation Manual, check the correctness and quality of assembly of the Device and its components, ensuring the operation in accordance with the established requirements.

Do not trust the production of these works to random people, avoid independent unskilled actions!

Remember that if you violate the rules of work, you may lose the right to free warranty repair!

ATTENTION!

The personnel who have studied this Manual and have experience in operating medical equipment are allowed to service the Device.

Please send your comments on the improvement of the design and usability to the address of the manufacturer or supplier.

CONTENT

1	General information	5
2	Purpose and scope.....	6
3	Composition and description of the Device ..	6
4	Basic security measures	18
5	Transportation	21
6	Placement requirements	22
7	Assembly and preparation for operation	23
8	Instructions for use of the Device	27
9	Maintenance, calibration and repair	40
10	Typical malfunctions and methods of their elimination.....	64
11	Disposal	66
12	Warranty	67
	Appendix «A»	68
	Appendix «Б».....	74
	Appendix «В».....	75
	Appendix «Г».....	104

1 GENERAL INFORMATION

1.1 This Operation Manual (OM) contains a technical description of the Automated device "Q-PROCESSOR" for "PROTIA" panels processing for in vitro diagnostics (hereinafter referred to as the Device), rules and instructions for its safe operation and other information that is necessary know users and personnel performing transport, adjustment, maintenance and repair.

1.2 This Manual applies to all modifications of the Device.

If necessary, a Supplement to the Operation Manual is issued for each modification, which is included in the list of operational documents.

1.3 The purpose of this Manual is to provide all the information necessary for the transportation, commissioning, operation and routine maintenance of the Device.

1.4 This Manual has been executed in accordance with GOST R 2.610-2019, GOST R 2.601-2019, GOST IEC 82079-1-2014, GOST R 56018-2014 / ISO / IEC Guide 37: 2012, GOST R 54329-2011, GOST R 2.105-2019, Order of the Ministry of Health of the Russian Federation of January 19, 2017 No. 11n "On approval of requirements for the content of technical and operational documentation of the manufacturer (manufacturer) of a medical device" and the Decision of the Council of the Eurasian Economic Commission dated November 10, 2017 No. 106 "On Requirements for the Implementation, Maintenance and Evaluation of the Quality Management System of Medical Devices, Depending on the Potential Risk of Their Use."

1.5 Terms and definitions - according to GOST R IEC / TO 60788-2009, GOST 31508-2012, GOST R 52002-2003, GOST 18322-2016, GOST 27.002-2015, GOST 25866-83, GOST R 55991.7-2014, GOST R 50267.0-92 and GOST R 2.610-2019.

1.6 This Manual should be considered as an integral part of the Device and should be always available and at the disposal of the installer, user (s) and operator (s) for repair and maintenance.

ATTENTION!

The Analyzer is serviced by personnel who have studied this Manual and have experience in operating medical equipment.

2 PURPOSE AND SCOPE

2.1 The device is intended for processing of PROTIA panels by the immunoblot method for subsequent quantitative determination of the concentration of allergen-specific IgE antibodies in human serum or plasma by photometry using the Device in medical institutions, as well as in medical and research laboratories for the diagnosis of allergies.

Functional purpose:

Diagnostic aid. The use of the product has no population and demographic aspects.

The intended users of the device are a clinical laboratory physician, a medical laboratory assistant, and other medical professionals.


2.2 The device is a non-invasive medical device for long-term use, full factory readiness, which includes all the necessary to ensure reliable operation of the device, monitoring and control.

2.3 The design of the Device complies with the design documentation and meets the specified requirements according to its intended purpose.

By type of placement the device is classified as desktop, closed.

2.4 Operating mode - automatic, repeated-continuous.

2.5 Device according to the classification of GOST 31508-2012, "Rules for the classification of medical devices depending on the potential risk of use" (approved by the Decision of the Eurasian Economic Commission dated December 22, 2015 No. 173) and Order of the Ministry of Health of the Russian Federation No. 4n dated June 06, 2012 "On approval of the nomenclature classification of medical products" (as amended on 25.09.2014) refers to medical devices of potential risk degree 2a, subgroups 5.1, and according to the possible consequences of failure - to products of class B GOST R 50444-2020.

3 COMPOSITION AND DESCRIPTION OF THE DEVICE

3.1 The main functional parts of the Device are shown in Figure 1.

The device works as part of a complex with an external monitor, keyboard, mouse.


3.2 Package Contents Devices: main unit - 1 pc., MC-4 / S power cord - 1 pc., power extension cord - 1 pc., 60 ml bottle - 2 pcs., 125 ml bottle - 2 pcs., DI water supply tube - 1 pc., container for collecting waste water - 1 pc., tube for ensuring waste water drain - 1 pc., tray adapter for Q64 panels - 1 pc., tray adapter for 96M panels - 1 pc., instruction manual - 1 pc.

3.3 The device provides laboratory research based on the immune responses of an antigen with an antibody.

The material is plasma or blood serum.

The method is highly sensitive and specific.

3.4 Operational and technical characteristics of the Device are shown in Table 3.1.

Picture 1

1) Name of the part: Power cord

Photo:

Raw materials: PVC, copper, steel.

Model name and part (catalog) number: Power cord MC-4/S.

Dimensions: 200 cm.

Weight: 190 g.

Quantity (PCS) : 1.

Manufacturer name : Elcom International Pvt Ltd

Address: 20, Ground Floor, Industrial Estate, 408, SVS Rd, Prabhadevi, Mumbai, Maharashtra 400025, India.

2) Name of the part: Power strip

Photo:

Raw materials: PVC, ABC-plastic, copper, steel.

Model name and part (catalog) number: F9E40.

Dimensions: Length : 260 mm X Width 50 mm X Height 35 mm.

Weight: 400 g.

Quantity (PCS) : 1.

Manufacturer name : Elcom International Pvt Ltd

Address: 20, Ground Floor, Industrial Estate, 408, SVS Rd, Prabhadevi, Mumbai, Maharashtra 400025, India.

3) Name of the part: Bottle with a capacity of 60 ml

Photo:

Raw materials: HDPE

Model name and part (catalog) number: Narrow mouth bottle 60 mL & 580110

Dimensions: ID neck : 13.5 mm, OD -bottle : 39.8 mm, Height without closure : 82 mm

Weight: 14.5 g


Quantity (PCS): 2

Manufacturer name: Tarsons

Address: 902, Martin Burn Business Park, BP-3, Saltlake, Sector-V, Kolkata-700091, India

4) Name of the part: Bottle with a capacity of 125 ml

Photo:

Raw materials: HDPE

Model name and part (catalog) number: Narrow mouth bottle 125 mL & 583120

Dimensions: ID neck : 17.7 mm, OD -bottle : 51 mm, Height without closure : 98.2 mm

Weight: 22 g


Quantity (PCS): 2

Manufacturer name: Tarsons

Address: 902, Martin Burn Business Park, BP-3, Saltlake, Sector-V, Kolkata-700091, India.

5) Name of the part: Tube for supplying DI- water

Photo:

Tube.

Raw materials: Tygon S3 E-3603 NSF 51

Model name and part (catalog) number: Tygon - ACF00004

Dimensions: ID : 2.4 mm x OD 4 mm

Weight: 4.5 g

Quantity (PCS): 1

Manufacturer name: Sanit Gobain

Address: 2664 Gilchrist Rd, Akron, OH 44305, USA

Connector.

Raw materials: Nylon

Model name and part (catalog) number: Luer Fitting, male Luer X 1/8 HB, White Nylon & LM4130

Dimensions : ID 2.3 mm x OD 4.8 mm

Weight: 0.5 g

Quantity (PCS): 1

Manufacturer name: Colder products company

Address: NORTH AMERICA, Roseville - 2820 Cleveland Ave. N., Roseville, MN 55113

6) Name of the part: Tube for the drain of wastewater

Photo:

Tube.

Raw materials: PVC NSF 51

Model name and part (catalog) number : PVC & K010-0204

Dimensions: ID : 3.2 x OD 6.4 mm

Weight: 47 g

Quantity (PCS): 1

Manufacturer name: Kuri tec

Address: 360 E. STATE PARKWAY SCHAUMBURG, IL 60173-5335 USA

Connector.

Raw materials: Nylon

Model name and part (catalog) number: Luer Fitting, male Luer X 1/8 HB, White Nylon & LM4130

Dimensions: ID 2.3 mm x OD 4.8 mm

Weight: 0.5 g


Quantity (PCS): 1

Manufacturer name: Colder products company

Address: NORTH AMERICA, Roseville - 2820 Cleveland Ave. N., Roseville, MN 55113

7) Name of the part: Tray adapter for Q64 panels

Photo:

Raw materials: Aluminium

Model name and part (catalog) number: Q 64 Tray Adapter & CMM0034

Dimensions: Width: 122 mm x Height: 90 mm

Weight: 64 g


Quantity (PCS): 01 number

Manufacturer name: CPC diagnostics Pvt Ltd.

Address: New No.: 324, Old No.: 237/1, Silver Jubilee Building, Khadi Gramodyog Complex, Avvai Shanmugam Salai, Goplalpuram, Chennai - 600086, India.

8) Name of the part: Tray adapter for 96M panels

Photo :

Raw materials: Aluminium

Model name and part (catalog) number : 96M Tray Adapter & CMM0035

Dimensions: Width: 122 mm x Height: 90 mm

Weight: 70.5 g

Quantity (PCS): 1

Manufacturer name : CPC DIAGNOSTICS PVT. LTD

Address: #9, V Floor, Gokul Tower, No. 9&10, C.P.Ramaswamy Road, Alwarpet, Chennai – 600018, India

9) Name of the part: Container for collecting waste water

Photo:

Raw materials: HDPE

Model name and part (catalog) number: Water tank 1 L “QC-18-01”

Dimensions: ID neck : 33 mm, outer sizes: (L) 100mm x (W) 64mm x (H) 230mm

Weight: 80 g

Quantity (PCS): 1

Manufacturer name : CPC DIAGNOSTICS PVT. LTD

Address: #9, V Floor, Gokul Tower, No. 9&10, C.P.Ramaswamy Road, Alwarpet, Chennai – 600018, India

Table 3.1

Parameter name	Value, norm
1	2
Appearance, installation quality	in accordance with the wiring diagram and reference samples; defects are not allowed
Rated supply voltage, V	220 ± 10%
Current type	variable single phase
Power supply frequency, Hz	50±5 %
Maximum current consumption, A, no more	1
Power consumption, VA, no more	200
Minimum sample volume for research, µl	50 or 100
Weight, kg, no more	30
Overall dimensions (length, width, height), mm	510 x 383 x 473 ± 5%
Climatic version	UHL 4.2
Group depending on mechanical stress	2
Sound power level, dBA, not more	60
Time of continuous work, h, not more	6 (followed by a break of at least 0.5 h)
Number of simultaneously processed panels	10
Interfaces for communication with a mouse, keyboard, monitor	USB 2.0, RS232
Time from loading the sample into the Device until the first panel is ready for scanning, min. no more	150
Performance (the number of panels ready to scan, obtained in one hour of the Device operation)	4
Heating time during incubation to a temperature of 37 °C after switching on, min., No more	45
Requirements for constructive compatibility with the computer	according to GOST 22315-77 and GOST 22316-77
Minimum monitor requirements	VGA 1024 × 768, 16M colors

i Note:

1 * *Permissible dimensional deviations - according to working drawings, GOST 24643-81, GOST 25346-2013 and GOST 25347-2013 (not worse than grade 14).*

2 *The manufacturer reserves the right to change and clarify the characteristics given in table 3.1 in accordance with the requirements of the design documentation for a specific modification of the Device.*

3.5 The device is equipped with software in accordance with GOST R IEC 62304-2013 ensuring its operation in the stated modes, compatible with operation system *Windows10*.

Software class according to GOST R IEC 62304-2013 - A.

The software performs the functions:

- control of the set operating modes;
- display of status and parameters of work;
- execution of the specified algorithms;
- diagnostics of equipment operability with the issuance of messages;
- operational management;
- registration of emergency events;
- formation of a database;
- generating reports;
- changes in control parameters, equipment configuration.

Pop-up windows, trends, graphs of changes in the measured technological parameters are used to display information.

3.6 It is possible to integrate the Device into a laboratory information system (LIS).

3.7 Device operating conditions

3.7.1 The device is suitable for operation in a UHL climate, location category 4.2 in accordance with GOST 15150-69

Permissible ambient temperature during operation: from plus 18 to plus 32 °C, relative humidity: not lower than 20% and not higher than 80% at plus 25 °C, atmospheric pressure: from 84.0 to 106.7 kPa.

3.7.2 Altitude above sea level: no more than 2,000 m.

3.7.3 Type of atmosphere - I according to GOST 15150-69.

Characteristics of the environment and aerosols contained in it: non-toxic, non-explosive according to GOST 30852.0-2002, GOST 30852.19-2002, 4-th hazard class according to GOST 12.1.007-76.

3.7.4 The device is resistant to mechanical factors in product group 2 in accordance with GOST R 50444-2020.

3.7.5 The device retains its characteristics when exposed to constant magnetic fields or alternating fields of mains frequency with an intensity of up to 400 A / m.

3.7.6 The device is allowed to be used only with chemical reagents, solutions, compatible with its design.

3.7.7 The permissible angle of inclination to any side during use is no more than 5°.

3.7.8 The device, when transported in the manufacturer's packaging, retains its properties after exposure to transport shock with a vibration frequency of 2 ... 3 Hz, an acceleration of 30 m*s⁻².

ATTENTION! The use of the Device in other conditions is considered as its use for other purposes.

3.8 Ergonomic and aesthetic characteristics of device generally complies with GOST 30.001-83, GOST 20.39.108-85, GOST R 56274-2014 and GOST R EN 614-1-2003.

Switches and switches meet the standards of GOST 23000-78.

3.9 Controls The devices are provided with inscriptions (symbols) in accordance with GOST 24505-80, GOST R IEC 60073-2000 and GOST R IEC 878-95, indicating region management.

Coloring and selection of signal colors are made taking into account the instructions of GOST 12.4.026-2015 and GOST R 50267.0-92.

3.10 The device has an indication

- power connections;
- the occurrence of a malfunction;
- data reception;
- readiness to perform operations;
- the end of operations.

The equipment contains instructions in Russian regarding:

- the main connections or related aspects;
- indication of the position of the main switch "on";
- the presence of high voltage;
- the data of rated current and characteristics of fuses;
- the conditions for opening the case.

3.11 Electrical equipment is installed in accordance with the requirements of the design documentation and the electrical schematic diagram according to RDT 25 106

3.12 The design of the Device is made controllable and maintainable in accordance with GOST 23660-79, GOST R 27.605-2013, GOST 26656-85, accessible for inspection and removal of components to be re-

placed, and (as necessary) their sanitization.

Repairs are carried out in specialized organizations.

3.13 The device has been designed so that its electronic part is maximally adapted to the use of automatic assembly.

The device contains a minimum number of elements for setting the electronic part, but not to the detriment of its technical parameters.

3.14 Surfaces Devices made of non-corrosion-resistant materials are coated in accordance with GOST 9.032-74, GOST 9.031-74, GOST 9.301-86 (class not lower than III) and GOST 9.104-2018 (group of conditions 4). Metallic and non-metallic inorganic coatings comply with GOST 9.301-86 and GOST 9.303-84.

3.15 Indicators of device reliability.

3.15.1 Mean time between failures of the Device is at least 800 hours.

The failure of the Device is a violation of its operational state associated with the failure of any component, which entailed deviation of the parameters beyond the declared limits, if at the same time, to restore the operational state, it is necessary to replace or repair any component.

3.15.2 The average time to restore the working condition when replacing a defective part (element) with a serviceable one should be no more than 2 hours.

Duration of daily maintenance - no more than 1 hour, monthly - no more than 5 hours, annual - no more than 10 hours.

3.15.3 Average service life of the Device is at least 5 years.

The criterion of the limiting state is considered to be technical and economic inexpediency restoration of characteristics.

3.15.4 Coefficient of technical use of the Device - not less than 0.90.

The degree of readiness is not less than 0.99 in accordance with GOST 27.002-2015.

3.15.5 The components, devices, apparatuses and devices used in the design have at the time of installation in the Device a residual service life not less than its warranty period.

3.15.6 The Device Software has:

- noise immunity, that is, it excludes the possibility of looping and getting into deadlock situations, the ability to function correctly in the event of failures, failures of some computing facilities, personnel errors, and also provides an automatic restart when the power supply is restored after a shutdown without giving false signals and control actions;

- speed, that is, the minimum restart time, reactions to external events, minimum loss of time for data protection;

- adaptability and development prospects, i.e., the simplicity of adapting programs to changes or extensions of the user's tasks without impairing other indicators.

3.16 The design of the device is suitable for disinfection, ensures the drainage of washing liquids; the number of grooves and protrusions that facilitate the accumulation of dirt or body fluids is minimal.

The construction of the Device does not contain any waterproof places, partitions, sharp narrowing of the cross-section and other inaccessible places.

3.17 Patent clearance for manufacturing methods and design solutions is ensured in relation to the Russian Federation and countries where the device can be delivered, as well as the transfer of technical, information and other documentation.

Joint intellectual property with commercial potential obtained in the course of work must be protected in accordance with Russian and international legislation.

3.18 The body of the Device has a smooth, uniform and uniform outer surface, without burrs, cracks, stains, dirt, scratches, traces of corrosion, cavities, dents and chips, in accordance with approved reference standards.

3.19 The materials used are approved for use in medical equipment by the bodies and institutions of Roszdravnadzor.

The materials are compatible with biological fluids.

3.20 By agreement with the customer, it is allowed to manufacture the Device with additional technical requirements for the composition, configuration and operating parameters, which should be indicated in the order.

3.21. Storage and operation The devices are carried out in closed, well-ventilated rooms in conditions that exclude moisture condensation, exposure to direct sunlight and aggressive environments (dust, acid and alkali vapours, as well as gases that cause corrosion).

3.22. The room should be free from aggressive environments, the mass concentration of dust in the air should not exceed 0.75 mg / m², the electric component of the electromagnetic field should not exceed 3 V / m in the frequency range from 150 kHz to 300 MHz.

3.23. The device complies with the requirements of GOST R IEC 61326-2-6-2014 for immunity and electromagnetic emission given in this standard.

- The device is designed and tested in accordance with the requirements of CISPR 11 (class A). In residential areas, it can cause radio interference, in which case the user should take measures to reduce the level of interference if the electromagnetic emission complies with class A.
- The user should assess the electromagnetic environment before using the device.
- Do not use the Device near sources of strong electromagnetic radiation (for example, unshielded intentional sources of radio frequency radiation), as they may interfere with its normal operation.

- If it is necessary to operate the Device near powerful sources of electromagnetic emission, or with equipment with extremely low noise immunity, take the following possible measures:
- time protection;
- distance protection;
- reduction of radiation parameters in the radiation source itself;
- rational placement of equipment in the laboratory with an account of the individual characteristics of each piece of equipment in terms of noise immunity and electromagnetic emission;
- rational modes of operation of installations and work of laboratory personnel;
- allocation of radiation zones.

**ATTENTION!**

It is not allowed to remove identification plates from the Device: they bear the serial number and other useful information.

Changes and alterations to the Device by the user, not agreed with the manufacturer, are not permitted.

4 BASIC SECURITY MEASURES

4.1 The device meets the safety requirements when used for the purposes specified in this manual.

The design ensures the reliability and safety of its operation during the established service life and provides for the possibility of maintenance, sanitary cleaning, repair and operational control.

4.2 During preparation for operation and operation of the Device, the measures of safety, labor protection and environmental protection specified in this Manual and the Operation Manual of "PROTIA" panels must be observed.

4.3 Safety instructions

4.3.1 The device complies with the structural and electrical safety standards

GOST IEC 60065-2013, GOST R ISO / TO 16142-2008, GOST R 50267.0-92 (according to class I, type B), GOST R 52770-2016, GOST IEC 61010-1-2014, GOST IEC 61010-2-081-2013, GOST IEC 61010-2-101-2013, GOST IEC 60950-1-2014.

4.3.2 Power cord MC-4 / S.

Composition: PVC, copper, steel.

Length: 200cm.

Weight: 190g.

Manufacturer's name: Elcom International Pvt., Ltd., India.

4.3.3 The temperature of the parts of the Device, accessible for touching without using a tool, during normal operation and at an ambient temperature from plus 18 to plus 32 ° C does not exceed plus 41 ° C in accordance with GOST R 50444-2020.

4.3.4 Sanitary and hygienic requirements for the Device are provided in accordance with the "Uniform Sanitary and Epidemiological and Hygienic Requirements for Goods Subject to Sanitary and Epidemiological Surveillance (Control)" dated May 28, 2010 No. 299, chapter II, section 18. The materials of the device do not provide irritating, toxic and sensitizing effects; no potentially hazardous radiation.

4.3.5 The design of the Device eliminates the possibility of contact with waste and the possibility of infection of the operator, and also prevents biological contamination.

Note - *During use, there is no contact of reagents and biological material with human skin or other organs, so there is no need to conduct toxicological and sanitary tests.*

4.3.5 The device meets the requirements of electromagnetic compatibility in accordance with GOST R IEC 61326-1-2014, GOST R IEC 61326-2-6-2014 and GOST R 51318.11-2006 (class A, group 1), as well as GOST IEC 61000-3-2-2017 (class A) and GOST IEC 61000-3-3-2015.

4.3.6 The design of the device components excludes the possibility of cuts and injuries during its use and repair.

The Device has no accessible metal parts galvanically connected to the power supply.

4.3.7 The design of control electronic devices excludes:

- not envisaged change of operating modes;
- the presence of errors during work;
- the presence of systematic failures.

Electrical equipment of devices are short-circuit and no-load resistant.

The device has 2 fuses of 2 amperes each.

4.3.8 Requirements for electrical safety (resistance, insulation strength, grounding and impedance) are provided in accordance with GOST R 50444-2020.

Permissible supply voltage deviations: ± 22 V in accordance with GOST R 50444-2020.

4.3.9 The equivalent sound level during operation of the Device does not exceed 60 dB in accordance with GOST 12.1.003-2014. Calculation methods - according to GOST 30530-97.

The vibration level during the operation of the Device does not exceed the standards of GOST 12.1.012-2004.

4.3.10 The degree of protection of electrical equipment is IP21 in accordance with GOST 14254-2015.

4.3.11 The Device software is protected against unauthorized configuration by third parties (password protection, encryption, access control) in accordance with GOST R 50739-95 and GOST R 52069.0-2013, as well as protection against incorrect use of functions and the input of invalid data.

The occurrence of a situation in which it is possible to access the functions of the Device without authorization has been excluded.

4.3.12 Requirements for fire safety The devices are in accordance with GOST IEC 60065-2013 and "Technical Regulations on Fire Safety Requirements" TR RF 005/2008 (Federal Law of the Russian Federation No. 123-FZ dated 22.07.2008).

The probability of a fire in (from) the Device is more than 10^{-6} per one product per year.

4.3.13 The electrical circuit excludes the possibility of spontaneous switching on, off and changing the operating modes of the Device; complete or partial interruption of power supply and subsequent restoration does not lead to hazardous situations, including spontaneous start-up.

4.4 During the operation of the Device, the following must be provided:

- periodic control over it;
- timely implementation of the regulations for prevention and repair;
- compliance with the requirements for safe placement.

4.5 In addition to the requirements of this Manual, during the operation of the Device, it is necessary to comply with the requirements of scientific and technical documents of factories - manufacturers of its electrical equipment and devices.

4.6 Repair of the Device must be carried out by the manufacturer or an organization authorized by him. Self-elimination of malfunctions (except those regulated in this Manual) and performance of repair and adjustment work is not allowed. It is forbidden to modify the design of the Device outside the factory conditions.

4.7 It is prohibited to use the Device if it does not comply with the manufacturer's passports, as well as the requirements of regulatory and technical documentation.

4.8 It is forbidden to use panels other than "PROTIA".

WARNING!

In case of violation of the norms and rules of operation, the requirements of safety measures established in this Manual, even if the violation was single and related only to one of the established requirements, rules and regulations, the manufacturer and seller, regardless of the timing of purchase and duration of operation

of the Device bear any responsibility for its quality and technical condition, as well as for any consequences that occurred during installation and / or during preparation for operation and / or during operation of the Device, including those that caused damage to the health and life of people, damage environment and human habitat

IT IS FORBIDDEN to use the Device:

- upon detection of mechanical and other damage to components, electrical wiring, fuses and switching devices;
- in case of damage to structural elements of the equipment.

IT IS FORBIDDEN to perform adjustment, repair and similar work with the electrical equipment switched on.

ATTENTION! Failure to comply with these and other safety measures and precautions specified in this manual may endanger the life and health of the user, and cause accidents.

5 TRANSPORTATION

5.1 Requirements for transportation - according to GOST 23216-78 and GOST R 51908-2002.

Before transporting the Device, make sure that the package is complete and that there are no violations in the packaging.

5.2 Component parts of the Device must be protected from movement by gaskets and shock absorbers made of corrugated cardboard.

Changeable, spare parts and accessories packaged in polyethylene packages.

5.3 Transportation of devices must be carried out by carriers specialized in the transportation of the relevant goods and having the proper licenses (permits) and experience in transportation.

5.4 The Device is shipped in a shipping container by any covered vehicle with protection against the effects of atmospheric precipitation and fog, including the following:

- direct transportation by road for a distance of up to 1,000 km on roads with asphalt or concrete surfaces (roads of the first category) without speed limits or at a speed of up to 40 km / h for a distance of up to 250 km on cobblestone or unpaved roads (roads of the second and third categories);

- multimodal transportation by rail, air (in heated pressurized compartments), river modes of transport, in combination with each other, and by road, as well as transportation by sea.

5.5 Conditions of transportation. The devices in terms of the impact of climatic factors must comply with group 2 (C) in accordance with GOST 15150-69, in terms of the impact of mechanical factors - group G in accordance with GOST 23216-78.

Transportation temperature should be within the range from -10 °C to +40 °C, relative humidity: from 20 to 80%, atmospheric pressure: 86... 106 kPa.

5.6 Loading and unloading of the Device should be carried out in accordance with GOST 12.3.009-76.

The placement and fastening of containers in vehicles should ensure its stable position, exclude the possibility of displacement and impacts of them against each other, against the walls of the shipping container and vehicles.

5.7 It is forbidden to stack boxes during transportation or place foreign cargo on them.

6 PLACEMENT REQUIREMENTS

6.1 The placement of the Device must be carried out taking into account the technical data and safety parameters specified in the technical documentation.

Storage conditions - in group 1 (L) according to GOST 15150-69.

6.2 Storage of devices are carried out in closed, well-ventilated rooms in conditions excluding moisture condensation, exposure to direct sunlight and aggressive environments (dust, acid and alkali vapors, as well as corrosive gases).

Storage temperature should be within the range from plus 10 to plus 40 °C, relative humidity: from 20 to 90%, atmospheric pressure: 86... 106 kPa.

6.3 The distance from the packaging placed on the racks to any object (including the walls and floor of the store) should be at least 0.1 m, the distance to the nearest heating device should be at least 1.0 m.

6.4 The device, transportation, use and repair of which is not planned within 3 - 30 days, should be put into short-term storage, and with a duration of more than 30 days. - for long-term storage.

6.5 Information about the storage of the Devices, the operating organizations must record in the form (passport), which is part of the operating documents, which indicate the inventory number, completeness, start date and removal from storage.

6.6 Control of the technical condition and safety of the Devices should be carried out at least once a month for short-term storage and once every 3 months for long-term storage.

7 ASSEMBLY AND PREPARATION FOR OPERATION

7.1 After unpacking the Device must be placed on the desktop.

ATTENTION!

If the Device has already been in use, it is necessary to check that it has a decontamination declaration attached to it in accordance with Appendix «A».

NOTE - *It is recommended not to throw away an empty transport box that may be required for future transportation, storage, etc.*

7.2 Next, you should make sure that the device is complete.

The delivery should include:

- main unit – 1 pc.;
- power cord MC-4/S – 1 pc.;
- power strip – 1 pc.;
- bottle with a capacity of 60 ml – 2 pcs.;
- bottle with a capacity of 125 ml – 2 pcs.;
- a tube for supplying deionized water – 1 pc.;
- a container for collecting waste water – 1 pc.;
- a tube to ensure the drain of waste water – 1 pc.;
- tray adapter for Q64 panels – 1 pc.;
- tray adapter for 96M panels – 1 pc.;
- operating manual – 1 pc.

7.3 When assembling the Device in order to ensure safety measures, it is prohibited:

- get down to work without examining the operational documentation;
- start work without checking the technical condition of the elements of the Device, its electrical cables and control devices;
- use elements, components, tools and devices not provided for in the design.

7.4 The Device should be installed in winter only after keeping it unpacked in a closed heated room for at least 4 hours.

7.5 The device is wired according to the electrical schematic diagram.

The safety of the installation of electrical equipment and components must be ensured by meeting the requirements of GOST 12.3.019-80.

7.6 The following requirements are imposed on the placement of the Device:


- it is prohibited to install it at an angle;
- the mains socket used to connect the mains plug must be in an accessible place to be able to quickly disconnect the Device from the mains;
- the device must be located at a distance of at least 1 m from heating devices;
- it is prohibited to switch on the Device without protective grounding;
- it is not recommended to place or throw heavy objects on the power cable and damage it in any way;
- it is not recommended to place the Device in places:
 - with mechanical vibration;
 - with high humidity;
 - with direct exposure to sunlight;
 - with a very low or high temperature;
- do not use the Device in the immediate vicinity of highly flammable materials;
- if it is necessary to move the Device to another place, at least two people should do this.

7.7 The room should be free from aggressive environments, the mass concentration of dust in the air should not exceed 0.75 mg / m², the electric component of the electromagnetic field should not exceed 3 V / m in the frequency range from 150 kHz to 300 MHz.

7.8 After the Device is installed in the desired position, it must be connected to a 220 V AC power supply using the supplied cables.

(To ensure correct operation of the instrument, the worktable must be connected to the laboratory's grounding system.)

7.9 Connecting personal computer and pipes (Figure 2)


Figure 2

1 Connect the keyboard, mouse and monitor to the KEYBOARD, MOUSE, MONITOR connectors, respectively.

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3 Connect the pipe to ensure the waste water drain to the WASTE WATER OUTLET.

4 Connect the power cable to the POWER jack.

7.10 The device must be fully assembled, installed and connected to the mains in accordance with GOST 32144-2013 and GOST R 50571.28-2006 in the following order:

- check the integrity of electrical elements and installation of control devices, the absence of twists in the cables; mount electrical equipment;
- connect all units of electrical equipment to the ground loop;
- check the absence of violations of the insulation of the electrical wiring;
- connect the equipment to the AC mains; when connecting, check the compliance of the mains voltage and electrical equipment;
- install and configure software;
- check the equipment by sequential switching on in the setup mode.

7.11 Preparation of samples and reagents

1 Necessary to check the components of "PROTIA" before testing for the following indicators:

- product code and test type;
- product shelf life;
- number of panels of allergens;
- the volume of each reagent.

2 All reagents should be warmed to room temperature and mixed thoroughly before use. The packaging of allergenic panels may only be opened after equilibration to room temperature.

ATTENTION!

Air bubbles are not allowed!

3 Prepare 1X washing solution by diluting concentrated 20x washing solution by deionized water.

4 Prepare test samples (serum or plasma).

5 Remove blood cells or any solids, if any, by centrifugation prior to testing.

ATTENTION!

Hemolyzed or contaminated samples may lead to incorrect results.

6 If frozen samples are used for testing, they should be completely dissolved, mixed and centrifuged before use (avoiding repeated freezing and thawing cycles).

7.12 Verification of operating modes (characteristics) is carried out when monitoring the operation of the Device.

7.13 Minimum monitor requirements: VGA 1024 × 768, 16M colors.

ATTENTION!

Incorrect or carelessly performed assembly, adjustment, preparation of samples and reagents can lead to the complete destruction of the Device or to the fact that it becomes unsafe for the user!

8 INSTRUCTIONS FOR USE OF THE DEVICE

8.1 Before starting operation, control of the correct assembly of individual equipment and "PROTIA"

panels in accordance with their documentation must be ensured.

Make sure by visual inspection that the Device is in good working order condition.

8.2 Operation should be carried out in accordance with SanPiN 2.1.3678-20 and GOST ISO / IEC 17025-2019 in clean rooms in accordance with GOST R ISO 14644-4-2002, GOST R ISO 14644-5-2005 and GOST R 52539-2006.

When working with the Device, one should observe the instructions of GOST R 53079.4-2008, GOST 27025-86, GOST R 52905-2007, GOST 31891-2012 and "Rules for the device, safety precautions, industrial sanitation, anti-epidemic regime and personal hygiene when working in laboratories (departments, departments) of sanitary and epidemiological institutions of the system of the Ministry of Health of the USSR", Moscow, 1981

8.3 During laboratory work, use rubber gloves in accordance with GOST 20010-93, goggles and special laboratory clothing (coat, suit or overalls).

8.4 The device is suitable for work with panels "PROTIA Allergy-Q Atopic panel", "PROTIA Allergy-Q 64 Food panel", "PROTIA Allergy-Q 64 Inhalant panel", "PROTIA Allergy-Q 96M panel" and other panels "PROTIA" (RU No. RZN 2021/14464 dated 28.05.2021).

8.5 "PROTIA" panel reference is presented in Appendix «Б» of this Guide.

8.6 Control of biocontamination in the laboratory should be carried out in accordance with GOST ISO 14698-1-2005 and GOST ISO 14698-2-2005.

ATTENTION!

The device involves the use of chemicals and other hazardous (caustic, irritating and hazardous) chemicals that can be hazardous to health.

The device provides for the processing of potentially infectious samples (plasma or blood serum) that may cause infection / contamination!

8.7 Algorithm for working with the Device

1 Place the reagents in the appropriate positions (Figure 1).

Place the PROTIA reagents in the reagent tray in the order from 1 to 4 in accordance with the color marking of the slots in the tray and the labels on the reagent vials according to Figure 3:

- sample diluent solution: # 1 - blue label;
- antibody solution: # 2 - yellow label;
- enzyme solution: No. 3 - red label ;
- substrate solution: # 4 - brown label.

Figure 3

ATTENTION!

Air bubbles from reagents should be removed prior to testing as they may affect the accuracy of the test.

2 Pour 1X Washing Solution into a 60 ml wash solution bottle and place the tubing into the wash solution bottle as shown in Figure 4.


Figure 4

3 Pour deionized water (DI) into a 125 ml DI bottle and place the deionized water intake tube into the DI bottle as shown in Figure 5.


Figure 5

4 Place the tray adapter with the panels installed in it into the tray of the Product (Figure 6). In this case, the black tray adapter must be used for the "PROTIA Allergy-Q 96M panel", and the transparent adapter for the rest of the panels.


Figure 6

5 Turn on the Device using the POWER BUTTON on the rear panel (Figure 2).

6 Double click the "Proteome" icon on the PC screen (Figure 7).


Figure 7

7 Press the "START" button on the main screen of the program (Figure 8).

Figure 8

8 Select test items and the number of allergen panels (1-10).
Then press the "CYCLE" button (Figure 9).

ATTENTION!

Select "ALLERGY-Q64" for all PROTIA Allergy-Q panels, except for the "PROTIA Allergy -Q 96M" panel (select "ALLERGY-Q96" for this).


Figure 9

WARNING!

If no product is selected, an alarm message appears as shown in Figure 10.


Figure 10

9 Confirm the adapter and reagents, and click the "OK" button (Figure 11).

Figure 11

10 Double-click the panel column you want to analyze and enter information about the sample. Click the "OK" button to go to the next step (Figure 12).

Figure 12

11 After pressing the START button, a new window will appear with a message like "Sample IDs are missing for 00 samples! Are you sure you want to continue? "

Here you should confirm the number of samples and click the "Yes" button as shown in Figure 13.

Figure 13

12 Check the information according to the messages that appear: "Make sure that a sufficient amount of reagents is placed in the appropriate places" and "Make sure the panels are loaded in the correct positions", then press the "OK" button sequentially to proceed to the next step (Figure 14).


Figure 14

13 The test will start immediately (Figure 15).

Figure 15

14 After the sample diluent solution has been added to each panel, a pop-up window will appear and an audible signal (“ringing of a bell”) will appear indicating that the test specimens (s) should now be applied by the technician (Figure 16).

To turn off the sound signal, press the "MUTE" button.


Figure 16

Samples should be entered manually, then click the "OK" button in the pop-up window as shown in Figure 16.

Subsequent processes are automatic.

16 When all processes are completed, a beep will sound informing you of the end of the test. To turn off the signal, press the "MUTE" button.

17 At the end of the process, the test panels must be allowed to dry, after which the results can be evaluated in accordance with the instructions for use of "PROTIA" (Appendix «Б» of this Manual).

8.8 Rinsing in mode "Prime"

1 Flushing of the Device should be performed after the end of each work session. To start it after completing the test, press the "Yes" button when the question "Do you want to fill the pumps?" (Figure 17).

2 When the message “Please place DI water in the position of the detergent solution” appears, you must replace the bottle with the detergent solution with the bottle of deionized water and press the “OK” button (Figure 18).


Figure 17


Figure 18

ATTENTION!

A flush cycle should be included at the end of each work session and after maintenance work.

3 When Prime mode ends, press the EXIT button to go to the main screen (Figure 8).

8.9 Shut Down the Device

1 Press the "Exit" button to end the program on the main screen (Figure 8).

2 Shut down MS Windows on the monitor (Figure 7) and turn off the power using the toggle switch on the rear panel.

8.10 Used panels "PROTIA" cannot be reused.

8.11 During the entire period of operation, the equipment of the Device must undergo a periodic technical examination at least once every 12 months.

After repair (replacement) or installation, the Device undergoes a partial technical examination in order to establish that the replaced, newly installed or repaired equipment is in good condition, ensuring safe operation.

8.12 Device upgraded during operation, undergoes a technical examination in the scope of a periodic check of a set of documentation for modernization, visual and measuring control.

8.13 A device that has served its service life is subject to examination in order to determine the possibility of extending its safe operation, including:

- visual and measuring control;
- control check of work in all modes;
- determination of the condition with the identification of defects, malfunctions, wear and corrosion;
- testing of safety devices;
- testing of protective grounding (grounding), insulation resistance of electrical networks and electrical equipment, checking the operation of protection in networks with a dead-grounded neutral.

Based on the results of the survey and analysis of operating conditions, work is being carried out to determine the residual resource of the device and the possibility of extending its safe operation period.

8.15 The owner of the Device must ensure its maintenance in good condition and safe operation by arranging proper maintenance. For this purpose, personnel must be appointed who are responsible for organizing and carrying out maintenance and repair work.

The maintenance personnel must know the device and the operation rules of the Device, pass the safety knowledge test and be admitted to independent work in the prescribed manner.

8.16 During the entire period of operation, it is obligatory to keep records of all work performed, to timely fill in the relevant sections of the form (passport) of the device.

ATTENTION!

- **the power cable should be disconnected if the Device is not used for a long period of time;**
- **when disconnecting the power cable, hold on to the plug and not pull on the cable itself;**
- **under no circumstances should you disassemble the Device yourself!**

- it is necessary to avoid getting water or any other liquid inside the Device case.

8.17 After each working session, it is necessary to check the cleanliness of the work table and, if necessary, thoroughly clean it and the corresponding accessories using gauze soaked in 0.5% sodium hypochlorous acid solution.

8.18 Decommissioning of the Device:

1 Pull out the appliance plug.

2 Thoroughly clean and disinfect the device following the disinfection instructions in Section 9 of this Manual.

3 Attach a completed Decontamination Declaration dated and signed by a qualified technician.

4 Pack for planned transportation, for which:

- place the device in the original packing box;

- affix appropriate labels to the outside of the packaging box with sufficient warnings to indicate that

"The product is out of service";

"The product has been in contact with potentially infectious biological substances."

5 Manual transport of the Device without packing it into the packing box must be confined to the working area and care must be taken not to damage the device.

8.19 During any shipment, the device must be accompanied by a Decontamination Declaration (Appendix B), preferably affixed to the packaging box.

8.20 Contraindications to use The devices are not regulated and are determined by the characteristics of the PROTIA panels.

Application restrictions are also imposed by the characteristics of the panels.

8.21 Specifications Devices stable in declared by 3.7 conditions exploitation.

Information on the traceability of values specified for calibrators or control materials, which is provided by means of available reference methods (methods) of measurement and (or) standards, is determined from the results of clinical trials.

8.22 Conditions required for collection, processing and preparation of samples, data on the stability of the analyzed samples, including conditions and duration of storage, transportation conditions, as well as imposed restrictions are determined by the characteristics of the PROTIA panels.

8.23 Taking medical decisions based on the results of testing without prior consultation with the attending physician is unacceptable!

9 MAINTENANCE, CALIBRATION AND REPAIR

9.1 Maintenance of the Device as a whole consists in regular inspection of all working units for serviceability and integrity, special attention should be paid to finding faults in electrical connections, damage

to protective coatings, the integrity of the body, insulation of electrical wiring, grounding devices.

Maintenance of the Device is carried out only after it is de-energized.

9.2 It is necessary to regularly visually check the operability and condition of the components and accessories of the device:

- free movement of the probe along the working area without obstacles and extraneous noise;
- free operation of the panel tray swing system without obstacles and extraneous noise;
- correct operation of the dosing and aspiration circuits (if necessary, replace the connecting tubes);
- no visible damage to the insulation of the power cord.

9.3. When using the Device, the external surfaces should be cleaned regularly, but at least once a week.

Use a cloth dampened with water and a mild detergent to clean all external parts. Do not use solvents of any kind (for example, alcohol or gasoline).

The bottles used must be kept clean.

Disinfectants are used as needed and must be ensured before application that they do not contain alcohol or abrasives.

NOTE - *Do not wet electrical parts (connectors, etc.) during cleaning. If electrical parts become accidentally wet, dry them thoroughly before turning on the appliance again.*

9.4 The operating time before calibration is from 3 to 6 months, depending on the intensity and conditions of use, or is carried out immediately in case of replacement of spare parts and assemblies.

9.5 Device Calibration

9.5.1 Precautions for Calibration:

9.5.1.1 Calibration should be performed on a flat surface with sufficient space around the instrument.

9.5.1.2 Do not move the Device during calibration.

9.5.1.3 Do not subject the Device to shock during calibration.

9.5.2 Preparing the Device

9.5.2.1 Make sure each component is properly connected to the Q-PROCESSOR [Figure 19].


9.5.2.2 Connect keyboard, mouse and monitor.

9.5.2.3 Connect the DI / deionized water tubing to the DI water supply.

9.5.2.4 Connect the waste water tubing to the Waste water outlet.

9.5.2.5 Connect the power cord to the Power Cord Connector.

9.5.3 Executing the service form

9.5.3.1 Double click the “Proteome” icon on the PC screen [Figure 20].


Figure 20

9.5.3.2 After completion of the test, all buttons will be activated [Figure 21].


Figure 21

9.5.3.3 Press the "SERVICE" button on [Figure 21], the "ENTER PASSWORD / ENTERPASSWORD" window will open [Figure 22].


Figure 22

9.5.3.4 If you enter the password code "eqw" and press the "OK" button, the "Service Form" will be displayed [Figure 23].

ServiceForm

X

Xhom Sensor XRef Sensor Zhom Sensor ZRef Sensor

Chom Sensor CRef Sensor Dhom Sensor

Probe Crash Sensor Fluid Sense Rocking hom Sensor

 Check

Fill/Empty Water Fill/Empty Buffer

Movement

Calibration 

Figure 23

9.5.4 Calibration of Probe Position

9.5.4.1 Click the Calibration button in the service form [Figure 23], the Calibration Form [Figure 24] will be displayed.

CalibrationForm

Positions

 Start

<< >> CHECK

GET VALUE SAVE CANCEL

Syringe Volume

 Syringe

Pump Volume

 Pump

Nozzle

 Start

<< >> Get Value

Cancel Save

Positions area

Syringe Volume area

Pump Volume area

Nozzle area

Figure 24

9.5.4.2 If you click the drop-down menu in the Positions area, you can check the calibration lists [Figure 25].


Figure 25

9.5.4.3 Select one of the calibration lists and press the Start button [Figure 25].

Example: "Common: Carousel"

9.5.4.4 The probe will move to the last position of its use [Figure 26].

Figure 26

9.5.4.5 If you press the "CHECK" button, the probe will move downward [Figure 27].


Figure 27

If the probe is not located in the middle of the reagent bottle, adjust its position using the "<<" or ">>" button.

Note: If you need to use "<<" or ">>" several times, press the button every two seconds.

9.5.4.6 After adjusting the probe position, press the CHECK button again and confirm the location (in the center of the reagent bottle) [Figure 27].

9.5.4.7 Press the GETVALUE button and the transmitter will move to the home position and the location will be recorded [Figure 28].


Figure 28

9.5.4.8 If you press the SAVE button, the message Calibration Saved [Figure 29] will be displayed and you can complete the calibration by pressing the OK button.

Figure 29

9.5.4.9 If you want to confirm the position after calibration, select Common: Carousel from the drop-down menu in the Positions area. Move the probe using the Start / Start button and press the CHECK button.

9.5.4.10 If after calibration the position of the probe is normal, press the "CANCEL / CANCEL" button, and the probe will move to the initial position.

9.5.4.11 If you are repeating processes 4.2 ~ 4.10, you can calibrate all lists in the drop-down menu. Each item in the drop-down menu is shown in [Figure 30] ~ [Figure 39].


[Figure 30] Common: Pump Calibration


Attention! In the Common: Pump Calibration list, DO NOT PRESS the CHECK button.


[Figure 31] Common: Nozzle Home


[Figure 32] Allergy-Q 64: Strip1_Probe / Probe in panel # 1


[Figure 33] Allergy-Q 64: Strip1_Pump / Aspiration from panel # 1

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
[Figure 36] Allergy-Q 96: Strip1_Probe / Probe in panel # 1

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[Figure 39] Allergy-Q 96: Strip6_Pump / Aspiration from panel # 6

9.5.5 Calibrating the movement of the nozzle (nozzle)

9.5.5.1 If you click the Calibration button in the service form [Figure 23], the Calibration Form [Figure 40] will be displayed.

The screenshot shows the 'CalibrationForm' window with a close button (X) in the top right corner. It is divided into three main sections, each outlined in red:

- Positions area:** Contains a 'Positions' label, a dropdown menu, a 'Start' button, a '<<' button, a '>>' button, a 'CHECK' button, a 'GET VALUE' button, a 'SAVE' button, and a 'CANCEL' button.
- Syringe Volume area:** Contains a 'Syringe Volume' label, a dropdown menu, a 'Syringe' button, and a 'Pump' button.
- Nozzle area:** Contains a 'Nozzle' label, a dropdown menu, a 'Start' button, '<<' and '>>' buttons, a 'Get Value' button, a 'Cancel' button, and a 'Save' button.

Figure 40

9.5.5.2 Click the drop-down menu in the Nozzle area and the calibration lists will be displayed [Figure 41].

This screenshot shows the 'Nozzle' area with the dropdown menu open, displaying two options: 'Dispense' and 'Aspirate'. A red arrow points to the 'Start' button. A red circular stamp is visible over the 'Dispense' option. To the right of the buttons, there is a red rectangular stamp with the text 'ProteomeTech Inc.' and '2010'.

Figure 41

9.5.5.3 Select one of the calibration lists.

For example: "Application of the solution / Dispense"

9.5.5.4 If you click the Start button [Figure 41], the message shown below [Figure 42] will appear.

Make Sure Allergy Q-64 Panel Is place on the tray


Figure 42

9.5.5.5 Place the panel on the tray and click the "OK" button. The nozzle then moves to the "Allergy-Q64: Strip1_Pump" location and is directed downward. [Figure 43].

Note: The aspiration nozzle should be approximately 2mm above the top of the panel.


Figure 43

9.5.5.6 If the position of the suction nozzle (nozzle) is too high or low, adjust the position using the "<<" or ">>" button.

Note: If you need to use the "<<" or ">>" button several times, press the button every two seconds.

9.5.5.7 Move the nozzle to the desired location and press the "GET VALUE" button. The probe will move to the home position and the location will be recorded.

9.5.5.8 If you press the SAVE button, the message Calibration Saved [Figure 44] will be displayed and you can complete the calibration by pressing the OK button.


Figure 44

9.5.5.9 If you repeat steps 5.2 ~ 5.8, you can calibrate all the dropdown menu lists in the Nozzle area. If you select the Dispense menu, the aspiration nozzle (nozzle) should touch the panel and the spring should be compressed by about 1-2 mm [Figure 45].


Figure 45

9.5.6 Calibrating the Unit Pump by Volume

9.5.6.1 Fill the solution bottle with distilled water and place it in the buffer bottle tray [Figure 46].


Figure 46

9.5.6.2 Insert the peristaltic pump tubing into the installed bottle.

9.5.6.3 If you click on the Fill / Empty Buffer button in the Service form, a warning will appear below [Figure 47].


Figure 47

9.5.6.4 If you click the "OK" button on [Figure 47], the message below will appear [Figure 48].


Figure 48

9.5.6.5 If you press the "OK" button on [Figure 48], the nozzle (nozzle) will move to its original position and begin to suck in distilled water [Figure 49].

Figure 49

9.5.6.6 When the message below appears after running the peristaltic pump, click the "OK" button and close the window [Figure 50].


Figure 50

9.5.6.7 Repeat processes 6.3 ~ 6.6 until the peristaltic pump tube is filled with bubble-free distilled water.

9.5.6.8 If you click the Calibration button in the Service form, the Calibration Form is displayed. Click the Pump button in the Pump Volume area [Figure 51].

Pump Volume



Figure 51

9.5.6.9 Then the nozzle (injector) moves to the "Common: Pump Calibration" position [Figure 52] and below the message [Figure 53] will be displayed.

Figure 52



Figure 53

9.5.6.10 Place the panel under the nozzle and click the "OK" button [Figure 54]. The Buffer Wash Pump will deliver distilled water.

Picture 54

9.5.6.11 After dosing, the probe will move to its original position and a pop-up window will appear for entering the measured volume value [Figure 55].


Figure 55

9.5.6.12 Enter the measured volume (μL) and press the "OK" button. The pump calibration factor will be calculated automatically and a message will appear below if the result is within the acceptable range [Figure 56].


Figure 56

9.5.6.13 Finish the calibration by pressing the "OK" button.

9.5.6.14 If the volume is out of range, a technical assistance message will be displayed below [Figure 57].


Figure 57

9.5.6.15 Click the "OK" button and contact the manufacturer or supplier for technical assistance.

9.5.7 Calibrating the Unit Syringe by Volume

9.5.7.1 Fill the deionized water bottle with deionized water and place on the back of the Unit [Figure 58].

Figure 58

9.5.7.2 Insert the tubing into the deionized water bottle.

9.5.7.3 If you click the Calibration button in the Service form, a warning message appears [Figure 59].


Figure 59

9.5.7.4 If you click the "OK" button on [Figure 59], the following message appears [Figure 60].


Figure 60

9.5.7.5 If you press the "OK" button in [Figure 60], the syringe pump will start and begin to suck in distilled water.

9.5.7.6 When the message shown below appears after the syringe pump is running, click the OK button and close the window [Figure 61].


Figure 61

9.5.7.7 Repeat processes 7.3 ~ 7.6 until the syringe pump tube is filled with bubble-free distilled water.

9.5.7.8 If you click the Calibration button in the Service form, the Calibration Form is displayed. Click the drop-down menu in the Syringe Volume area to view the calibration lists [Figure 62].


Figure 62

9.5.7.9 Select the scope of calibration.

For example: "250 µl"

9.5.7.10 If you press the Syringe button, the probe will move to the Common: Carousel position [Figure 63] and the message below [Figure 64] will appear.

Figure 63


Figure 64

9.5.7.11 Place the tube in the carousel position and click the "OK" button [Figure 64]. Then the syringe pump will deliver distilled water.

9.5.7.12 After dosing, the probe will move to its original position and a pop-up window will appear for entering the measured volume of water [Figure 65].


Figure 65

9.5.7.13 Enter the measured volume (μL) and press the "OK" button. The pump calibration factor will be calculated automatically and a message will appear below if the result is within the acceptable range [Figure 66].


Figure 66

9.5.7.14 Finish the calibration by pressing the "OK" button.

9.5.7.15 If the volume is out of the allowable range, a technical assistance message will be shown below [Figure 67].


Figure 67

9.5.7.16 Click the "OK" button and contact the manufacturer or supplier for technical assistance.

9.6 Control of the technical condition is carried out as part of laboratory equipment within the time limits established in the medical and prophylactic institution

9.7 Parts that are subject to wear or are defective and need to be replaced are requested as indicated:

- customer order numbers;
- type and / or version of the Device;
- its serial number;
- component code numbers;
- component descriptions;
- quantity;
- the name and address of the medical institution, and the delivery address.

ATTENTION!

Replacement of parts with original spare parts guarantees an efficient and long service life of the Device.

9.8 Specially trained personnel who have been instructed in accordance with GOST 12.0.004-2015, tested for knowledge of the Device design and have electrical safety qualification at least III are allowed for maintenance and repair.

ATTENTION!

Users can only perform simple routine operations (see decontamination instructions) of maintenance. In the event of any malfunctions or breakdowns caused by the use of the Device, corrective maintenance must be performed by a qualified specialist!

The manufacturer accepts no responsibility for unprofessional repair attempts and the resulting consequences!

9.9 During maintenance and repair work, do not smoke and use open fire, leave contaminated rags and other cleaning materials near and on the equipment of the Device.

10 TYPICAL MALFUNCTIONS AND METHODS OF THEIR ELIMINATION

10.1 Before contacting the service center, it is recommended that you study this section for the subject of determining the possible cause of the problem. A few simple checks and adjustments can eliminate it and lead to normal operation.

10.2 Possible malfunctions and methods of their elimination are shown in Table 10.1.

ATTENTION!

If the above methods of solving the arisen malfunctions do not help, you must contact the service department!

Unqualified actions can lead to injury, damage to the Device and cancellation of the warranty!

Performing any operation explicitly prohibited in this Manual, as well as any adjustments, assembly actions not recommended or prohibited in this Manual, will void the warranty!

Table 10.1

Troubleshooting of Q-PROCESSOR

No.	Failure (phenomenon)	Probable cause	Solution
1	The operating system does not load	The last trip was wrong.	Press and hold the reset button (on the back of the instrument) for 5 seconds and then release.

2	There are no icons on the desktop. Only the taskbar line is visible	Display signals can be extended to the next screen.	Press and hold the Windows key + press and release the P key until the project screen appears. Then use the mouse to select the Duplicate icon.
3	Probe failure reported	The reagent bottle is not open.	Open the cap of the reagent bottle.
4	No liquid detected.	Not enough reagent loaded	Add the appropriate reagent.
5	After aspiration of the spent solution, its residues remained in the panel	Blockage in the drain pipe	If the drain tube is bent sharply, make it straight. There must be a gap between the level of the waste fluid and the tip of the drain tube.
6	Probe cannot capture and hold reagent	Deionized water less / empty	Top up with deionized water.
7	The cleaning solution is fed at an angle	Blockage in the nozzle of the detergent solution dispenser.	Remove dirt in the nozzle.
8	Communication error	Controller board not powered	If the internal power control switch is off, turn it on.
9	The sensor panel is not aligned correctly with the panel.	Wrong panel loaded on swinging platform, or wrong analysis selection	Check the selected analysis. Loaded the desired panel.
10	The probe moves to the lid of the reagent bottle, returns to its original position and reports "Probe crashed."	The probe has moved up and is stuck outside the probe collision sensor.	Push the probe slightly and release it.
11	OS loading takes a long time	Windows system files have been corrupted	Run Command Prompt in Administrator Mode and run the "SFC / scan now" command.

11 DISPOSAL

11.1 In the event of the occurrence of limiting states and a decision on the unsuitability of the Device for repair and further operation or the inexpediency of further operation, its equipment must be dismantled and disposed of.

11.2 Before disposal, individual devices must be sorted out in order to assess the possibility of further use outside the device.

11.3 After the end of the service life, if further operation is impossible, the components of the Device, after dismantling, must be used or disposed of in the prescribed manner in specialized organizations.

The body and other metal parts of the Device are recyclable as secondary resources of ferrous and non-ferrous metals.

11.4 At the end of its service life in medical institutions, the Device is disposed of in accordance with SanPiN 2.1.3684-21 (waste class B), SP 2.1.7.1386-03 and Federal Law "On the Sanitary and Epidemiological Welfare of the Population" No. M 52-FZ dated March 30 .1999

11.5 Reagents and test media should be removed (discarded).

Waste materials that remain after analysis and cleaning are hazardous, carrying the risk of infection / contamination, therefore, their handling and disposal must comply with the requirements of the current regulations.

11.6 The personnel carrying out the disposal must comply with the safety rules for carrying out a particular method of disposal.

12 WARRANTY

12.1 The manufacturer guarantees the compliance of the Device with the requirements of technical and design documentation, subject to the rules of installation, operation, transportation and storage.

12.2 The warranty period of operation is 12 months from the date of receipt by the consumer (customer). The guaranteed shelf life is 1.5 years from the date of manufacture.

12.3 During the warranty period, the manufacturer carries out warranty repairs of the Device or equipment that has failed.

After repair and / or warranty service, the Device will be returned to the owner in its original condition at the time of delivery.

12.4 The warranty period for components and spare parts is set equal to the warranty period for the operation of the Device and expires simultaneously with the expiration of its service life.

12.5 Complaints are submitted during the warranty period, provided that the consumer complies with the requirements of the instructions for using the Device, with the preparation of a complaint report.

12.6 The warranty does not apply to:

- for scratches, normal wear and tear, wilful damage or damage due to negligence or non-observance of the instructions for use, improper installation or use of inappropriate voltage;
- damage due to chemical or electrochemical reaction, rust, corrosion or exposure to water; for damage due to being in an unusual (inappropriate) environment;
- for accidental damage caused by foreign objects or substances, as well as due to cleaning, for repairs carried out by unauthorized service organizations or unauthorized service partners; or repair using non-original spare parts;
- on cables, fuses, clamps and electrical connectors;
- at self-repair experiments without the manufacturer's permission;
- in case of damage caused by force majeure (typhoon, fire, flood, hostilities, etc.);
- in cases where the Device is used by unauthorized persons.

12.7 Submission to the Service Center of requirements to eliminate the defects of the Device is possible only with the simultaneous presentation of a correctly completed warranty card. At the same time, it must legibly indicate information about the Device, the defects in which are to be eliminated (name, serial number), as well as information about the sale of the Device (date of transfer to the buyer, name and address of the seller), certified by the signature and seal (stamp) of the seller, as well as the buyer's signature.

12.8 The manufacturer reserves the right, without notice, to make minor changes (modifications) to the design of the Device that do not affect the performance of the product as a whole.

APPENDIX «A»

Decontamination of Device

A.1 Decontamination (disinfection)

The use of the Device involves the processing of potentially infectious samples (plasma or blood serum) and positive control samples that may cause infection.

WARNING!

For the above reason, all parts of the instrument must be considered potentially contaminated.

A.2 Decontamination procedure

For decontamination, use one of the following broad spectrum disinfectants (or an equivalent):

- Lyseptol (Liseptol)
 - Aseptisol (Aseptisol)
 - Germocid Plus (Germanicide Plus),
- or similar ones, approved for use in medical institutions.

ATTENTION!

Disinfectants may only be used on the surfaces of the Device!

The following steps should be followed during disinfection:

- 1 Put on disposable gloves, safety glasses and suitable clothing...
- 2 Prepare an autoclave bag for all disposable items that need to be disinfected and label it with autoclave tape...
- 3 Disconnect the plug to disconnect the appliance from the mains and prevent short circuits....
- 4 Remove all accessories. All accessories that come with the device must also be disinfected.
- 5 Gently spray the surface disinfectant. Alternatively, wipe the entire instrument with gauze or soft paper dampened with disinfectant solution.
- 6 After a contact period of at least 10 minutes, repeat the previous steps (1-5) again.
- 7 After 5 hours of contact, clean the surfaces of the instrument with a mild detergent or deionized water to remove any disinfectant residue or dirt.
- 8 Thoroughly dry the Device.

- 9 Pack the Device and all accessories in the original packaging.
- 10 Disinfect your hands and wash them with a mild detergent.
- 11 Fill out the decontamination protocol and attach it to the device

A.3 Decontamination Declaration

The declaration must be completed and attached to each instrument before storage or before being delivered to a service center for maintenance or repair.

The declaration must be affixed to the packaging of the device in a conspicuous place.

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OVERVIEW: The following is the procedure for decontaminating spare parts in Q-PROCESSOR.

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Please find below the service policy for CPC Diagnostics Pvt Ltd. instruments that are potentially biologically contaminated and the decontamination procedure to be followed before returning instruments or parts.

All parts that are used in the sampling pathway may come in contact with human serum or other body fluids and must be decontaminated.

I. Introduction:

It is our policy to assume that any spare parts that come into contact with serum or other body fluids are potentially biologically contaminated. Field personnel should apply Universal Precautions whenever pos-

sible exposure to blood borne pathogens. Universal Precautions is an approach to infection control based on the concept that all human blood and certain body fluids are treated as if they were known to be infectious to HIV, HBV and other blood borne pathogens. Therefore, we assume that some spare parts are potentially biologically contaminated if the customer did not use samples of human serum or body fluids.

We recommend disinfectants: isopropyl alcohol 70% and bleach 5%.

The person responsible for the return of the parts listed below must carry out the decontamination procedure and complete the decontamination certificate.

II. Disinfection procedure:

1) Equipment:


If there is a possibility of biological contamination, the following minimum personal protective equipment must be used:

Gloves: must be worn when working on the surface of the instrument or anywhere within the sample path.


Safety glasses: When servicing tools that may spray liquid, such as liquid handling systems, wear safety glasses when working within the liquid path.

2) Spare parts requiring decontamination:

The following spare parts must be decontaminated before being returned:

Description	Product code	De-activation liquid	Photo
Waste water drain tube	SAG0153	Bleach 5%	
Tray adapter for 96M panels.	CMM0035	Isopropyl alcohol 70%	

Panel tray adapter Q64	CMM0034	Isopropyl alcohol 70%	
Waste flow control valve	CVA0010	Отбеливатель 5%	
Probe	CMR0028	Отбеливатель 5%	
Assay aspiration nozzle	SAG0247	Отбеливатель 5%	

3) Decontamination procedure:

- Wear gloves, a lab coat and safety glasses when performing the decontamination procedure.
- Apply disinfectant to the part to be disinfected.

Note: Be careful not to apply excessive amounts of disinfectant to the anodized parts, as this can remove small amounts of dye from the anodized parts.

- Leave the disinfectant on the surface for 10 minutes.
- Remove the disinfectant by wiping it with a damp cloth.

APPENDIX «Б»

(reference)

Packaging Details

The device is packed in a single-turn shipping container, which is a cardboard box 80x60x84 cm made of 7-layer cardboard, with subsequent installation of the shipping container on a pallet with interception

tape and protective film (Figure B.1). The device, together with the lower and upper cushioning pads made of polyethylene foam, is placed in a shipping container. The box fits inside the groove in the upper polyethylene foam pad.

Figure B.1

APPENDIX «B»

Instructions for use

Panels "PROTIA"

PROTIA *Allergy-Q* Atopic panel

In vitro diagnostic kit for the quantitative determination of allergen-specific IgE antibodies in human serum or plasma using the immunoblotting method

[Appointment]

The reagent kit is intended for the quantitative determination of allergen-specific IgE antibodies in human serum or plasma using the immunoblotting method.

[Functional purpose]

Diagnostic aid. The use of the product has no population and demographic aspects.

[Application area]

Clinical laboratory diagnostics.

Medical device for in vitro diagnostic use only.

[Potential consumers of the product]

Medical facilities of any profile.

[Professional level of potential users]

Clinical laboratory diagnostics doctor, medical laboratory technician.

[Possible side effects]

Side effects have not been identified.

[Indications for use]

Allergy diagnostics.

[Contraindications for use]

Severe hemophilia.

The product contains no materials that come into direct or indirect contact with the patient's body and the personnel using the product when meeting the requirements of the operational documentation (instructions for use).

The product contains no medicines and pharmaceutical substances.

The allergen panel contains components of animal origin:

- monoclonal antibodies against immunoglobulin E, Cat. No. 100105 (source of origin: mouse hybridomas, according to the documentation of the company Medix Biochemica (Finland));
- biotin-bovine serum albumin (source of origin: bovine serum albumin Cat. No. LY-0081 (bovine serum), according to the documentation of ProteomeTech Inc. (Republic of Korea)).

The sample diluent contains components of animal origin:

- gelatine Cat. No. G6650 (source of origin: bull flank - Type B, according to the documentation of SIGMA ALDRICH (USA)).

The antibody solution contains components of animal origin:

- biotin-conjugated mouse antibodies to human IgE, Cat. No. 100429/100106 (source of origin: mouse hybridomas, according to the documentation of Medix Biochemica (Finland));
- bovine serum albumin, Cat. No. LY-0081 (origin: bovine serum, according to the documentation of Boval Biosolutions LLC (USA)).

The enzyme solution contains components of animal origin:

- Bovine serum albumin, Cat. No. V5591 (source of origin: bovine serum), according to the documentation of the company Promega (USA));
- bovine serum albumin, Cat. No. LY-0081 (source of origin: bovine serum, according to the documentation of Boval Biosolutions LLC (USA)).

The components included in the kit are non-toxic and safe in the stated concentrations. Materials of animal origin are well fixed on the membranes of the test strip, placed in a plastic container with holes, which excludes contact with the human body when the medical device is used as intended.

Solutions, which include materials of animal origin, are well sealed with lids, when using solutions, a pi-

ette and gloves are used, which exclude contact with the human body when using a medical device for its intended purpose.

Thus, there are no components in the medical device that are toxic, aggressive and dangerous when using the medical device for its intended purpose and subject to basic safety measures.

Sampling, processing, storage and handling of these materials in accordance with the documentation of the supplier / manufacturer of these materials, as well as general safety rules for working with biological materials.

[Description of the target analyte, information about its scientific validity, principle of action]

Operating principle

PROTIA Allergy-Q is a multiplex diagnostic kit for allergy diagnostics based on the principle of enzyme immunoassay, which includes nitrocellulose membranes with different allergens adsorbed on regular interval lines, which allows you to test dozens of specific allergens in one test. PROTIA Allergy-Q uses a new parallel membrane technique and allows a wider variety of allergens to be assessed in a single test than other products using only one membrane. Allergen-specific IgE antibodies present in patient samples bind to antigens and are immobilized on the membrane even after the washing step.

Immobilized IgE antibodies bind to biotin antibodies to human IgE, and then biotin binds to alkaline phosphatase-conjugated lecnavidin. The color develops after the addition of the substrate in the last incubation with the enzyme, and the color intensity is analyzed with a color measuring device.

Description of the target analyte, scientific validity

Immunoblotting - a highly sensitive method for detecting proteins, based on a combination of electrophoresis and enzyme-linked immunosorbent assay (ELISA).

The pathogen antigens are separated by electrophoresis in polyacrylamide gel, then transferred from the gel to activated paper or nitrocellulose membrane and developed using ELISA. On these strips ("blots") the patient's serum is applied. Then, after incubation, the unbound antibodies of the patient are washed off and the serum against human immunoglobulins labeled with the enzyme is applied. The complex formed on the strip [antigen + patient's antibody + antibody against human Ig] is detected by adding a chromogenic substrate that changes color under the action of an enzyme.

Atopic allergy is a hypersensitive immunological condition mediated by immunoglobulin E (IgE) antibodies. Lymphocytic B cells, stimulated by a specific allergen, produce IgE antibodies to the antigen. IgE antibodies bind to receptors on the surface of mast cells or basophilic leukocytes through the Fc region. Subsequent binding of the allergen to cell-bound specific IgE causes cell degranulation and the release of vasoactive amines, resulting in smooth muscle contraction, itching, edema and trans-mucosal extracellular fluids. The most common clinical manifestations of this biological process are hay fever, asthma, dermatitis, urticaria, and anaphylactic shock.

[Kit Components]

Name of component	Characteristic	Quantity
1	2	
Allergen panel	Plastic cassette, with a plastic frame, containing a nitrocellulose membrane with adsorbed allergens (S-1, S-2, S-3, S-4, S-5 - standard lines; C - control line; from 1 to 44 - lines allergens)	20 panels

Sample diluent	A 20 ml polymer bottle with a blue screw cap containing 10 ml of blue-violet liquid	10ml x 1
Antibody solution	A 20 ml polymer bottle with a yellow screw cap, containing 10 ml of a yellowish or light yellow liquid	10ml x 1
Enzyme solution	A 20 ml polymer bottle with a red screw cap containing 10 ml of a colorless or yellowish liquid	10ml x 1
Substrate solution	A 20 ml polymer bottle with a black screw cap, containing 10 ml of yellow or yellowish liquid	10ml x 1
Washing solution 20x	A 20 ml polymer bottle with a white screw cap, containing 10 ml of concentrated colorless liquid	10ml x 1
Instructions for use		1

• Additional equipment required

1. Manual method: IMMUNOLOGICAL PHOTOMETRIC ANALYZER "Q-SMART" FOR IN VITRO DIAGNOSTICS ON "PROTIA" PANELS.

2. Semi-automatic method: "Q-PROCESSOR" AUTOMATED DEVICE FOR PROTIA PANEL PROCESSING FOR IN VITRO DIAGNOSTICS.

3. Automatic method: AUTOMATIC IMMUNOLOGICAL PHOTOMETRIC ANALYZER "Q-STATION ELITE" FOR IN VITRO DIAGNOSTICS ON PROTIA PANELS.

[Appearance]

1. Panel for allergens: a plastic panel in which a membrane is fixed in a plastikova frame.
2. Sample diluent: blue-violet liquid.
3. Antibody solution: yellow or light yellow liquid.
4. Enzyme solution: colorless or light yellow liquid.
5. Substrate solution: light yellow liquid.
6. Washing solution 20x: colorless liquid.
7. Instructions for use.

[Analysis methods]

1. Preparation of reagents and samples

(1) Preparation of the cleaning solution.

Dilute the detergent solution 20x with deionized water 20 times before use.

(2) Sample preparation

The test uses human serum or plasma (with anticoagulants sodium citrate, heparin, EDTA). Before testing, remove blood cells or any solids by centrifugation for 10-15 minutes at 2500-3000 rpm. Hemolyzed or contaminated samples may lead to incorrect results; therefore, it is not allowed to use samples with a high lipid content and (or) with signs of hemolysis, and (or) with visible microbial growth. Store serum or plasma samples at 2 ~ 8 °C when used for a short period (within 2 weeks), and at -15 °C or lower for longer storage. Thawed samples must be thoroughly mixed before testing. When working with thawed samples, it should be borne in mind that each thawing-freezing cycle leads to the destruction of detected antibodies, reducing the quality of the test material. In order to obtain reliable results of the study, it is recommended to use only a single freeze-thaw of the samples.

2. Analysis procedure

Manual method

1) Optimal conditions for the laboratory - conducting the experiment at $+ 25 \pm 2$ °C and at $50 \pm 20\%$ humidity.

2) All reagents should be brought to room temperature 30 minutes prior to use and mix thoroughly immediately before use. Open the packaging of the allergen panel after adjusting it to room temperature.

3) Fully moisten the test membrane by adding 300 µl of diluted detergent solution to the allergen panel, shake the panel for 5 minutes, and then remove any remaining detergent solution (100 rpm recommended).

4) Pipette 250 µL of Sample Diluent into the Allergen Plate.

5) Pipette 50 µL of patient sample into the panel and incubate with shaking at room temperature for 45 minutes.

6) Remove residual solutions from the panel and flush the membrane twice with diluted detergent solution. With each wash, add 300 µl of diluted wash solution, incubate with shaking for 5 minutes and drain. The solution should not remain in the panel.

7) Dispense 250 µl of antibody solution to the panel and incubate with shaking at room temperature for 30 minutes.

8) Flush the membrane according to method (6).

9) Pipette 250 µL of enzyme solution into the panel and incubate with shaking until for 30 minutes.

10) Flush the membrane according to method (6).

11) Dispense 250 µl of substrate solution and incubate with shaking at room temperature in a dark room for 20 minutes.

12) Remove the substrate solution from the panel and rinse it with 250 µl deionized water.

13) Leave the panel on the table to dry naturally (make sure the membrane completely dry).

14) Insert the panels into the IMMUNOLOGICAL PHOTOMETRIC "Q-SMART" ANALYZER FOR IN VITRO DIAGNOSTICS ON THE "PROTIA" PANELS and evaluate the results. For more details, see the Manual for IMMUNOLOGICAL PHOTOMETRIC ANALYZER "Q-SMART" FOR DIAGNOSTICS IN VITRO ON THE PANELS "PROTIA".

* to incubate and shake the panel, use any orbital shaker with a rotation speed (oscillation frequency) of 100 rpm. There are no restrictions on the choice of the device model.

Semi-automatic method

1) AUTOMATED DEVICE "Q-PROCESSOR" FOR PROCESSING PANELS "PROTIA" FOR DIAGNOSTICS IN VITRO is equipped with panels and reagents for allergens.

2) Refer to the operation manual AUTOMATED "Q-PROCESSOR" DEVICE FOR PROTIA PANEL PROCESSING FOR IN VITRO DIAGNOSTICS and let the reactions progress.

3) Insert the panels into the Q-SMART analyzer and evaluate the results. For more details, see the Manual for the IMMUNOLOGICAL PHOTOMETRIC ANALYZER "Q-SMART" FOR DIAGNOSTICS IN VITRO ON THE PANELS "PROTIA".

Automatic method

1) ANALYZER AUTOMATIC IMMUNOLOGICAL PHOTOMETRIC "Q-STATION ELITE" FOR DIAGNOSTICS IN VITRO ON PANEL "PROTIA" is equipped with panels and reagents for allergens.

2) Refer to the operation manual. AUTOMATIC IMMUNOLOGICAL PHOTOMETRIC ANALYZER "Q-STATION ELITE" FOR DIAGNOSTICS IN VITRO ON THE PANELS "PROTIA" and give the reactions progress to eat.

3. Evaluation and interpretation of results

1) The amount of allergen-specific IgE antibodies is determined in IU / ml using the Q-SMART Analyzer or the Q-STATION ELITE Analyzer, the class is determined from the table below.

2) Total IgE (tIgE) is expressed in IU / ml and can be quantified in the range of 1 to 2000 IU / ml.

Allergen-specific IgE		
ME / ml IgE	Class	Number
<0.15 ~ 0.34	0	Not found
0.35 ~ 0.69	1	Small
0.70 ~ 3.49	2	Moderate
3.50 ~ 17.49	3	Moderately large
17.50 ~ 49.99	4	Large
50.00 ~ 99.99	5	Very large
> 100	6	Extremely large

4. Quality control

The reference line (C) must be clearly visible. If the reference line is faintly visible, it cannot be read with optical measuring instruments. Retesting is recommended.

5. Characteristics

1) Analytical sensitivity (IgE detection threshold): 0.15 ME / ml.


2) Analytical specificity (interference and cross-reactions): None Interference and cross-reactivity of the kit with the following materials were found: bilirubin, albumin, immunoglobulins (IgA, IgM, IgG or IgD), 3 times higher than normal physiological levels.

3) Diagnostic sensitivity: 99.40% - 100% (95% confidence level).

4) Diagnostic specificity: 98.77% - 100% (with a confidence level of 95%).

5) Reproducibility: 100%.

6) Coefficient of Variation (CV): less than 8%

No.	Allergen	Code	Group
1	Total IgE	tIgE	IgE
2	Egg white	f1	Food
3	Egg yolk	f75	
4	Milk	f2	
5	Alpha lactoalbumin	f76	
6	Beta lactoglobulin	f77	
7	Casein	f78	
8	Peanut	f13	
9	Green pea	f12	
10	Soya beans	f14	
11	Walnut	f256	
12	Pork	f26	

13	Beef	f27	
14	Chicken meat	f83	
15	Cod	f3	
16	Mackerel	f206	
17	Crabs	f23	
18	Shrimp	f24	
19	Sea clam	f207	
20	Silkworm pupa	-	
21	Wheat	f4	
22	Baker's yeast	f45	
23	Rice	f9	
24	Potatoes	f35	
25	Peach	f95	
26	Apple	f49	
27	House dust	h1	Indoors
28	<i>D. pteronyssinus</i>	d1	
29	<i>D. farinae</i>	d2	
30	Cockroach	i6	
31	Cat epithelium and dandruff	e1	
32	Dog dandruff	e5	
33	<i>Candida albicans</i>	m5	Microorganisms
34	<i>Tricophyton rubrum</i>	m205	
35	<i>Aspergillus fumigatus</i>	m3	
36	<i>Penicillium notatum</i>	m1	
37	<i>Alternaria alternata</i>	m6	
38	<i>Staphylococcal enterotoxin B</i>	m81	Pollen
39	Birch-alder mixture	tx	
40	Oak	t7	
41	A mixture of cereal herbs	gx	
42	Ambrosia wormwood high	w1	
43	Common wormwood (Chernobyl)	w6	
44	Hop	w22	

[Precautions for Users]

- 1) For in vitro diagnostic use only.
- 2) PROTIA Alergy-Q can be used to aid in clinical diagnosis, the result is preliminary. The obtained analysis result should be correlated with the clinical picture. To clarify the diagnosis, a doctor's consultation is required.
- 3) False positive results can be obtained due to cross-reaction the activity of antigens tested with other antigens.
- 4) Do not smoke, eat or drink in areas where samples or tester are used new reagents.
- 5) Optimal conditions for a laboratory are to conduct an experiment at $+ 25 \pm 2$ °C and at $50 \pm 20\%$ humidity.
- 6) All specimens may include unknown infectious materials. When addressing. When handling samples, wear disposable gloves and wash your hands after testing.
- 7) To avoid accidents, do not place needles, knives, or any other object that could cause harm when handling specimens and reagents.
- 8) All patient samples and components used in the kit must be considered. Treated as sources of biohaz-

ards. They should be disposed of in accordance with the relevant guidelines.

9) Do not use kits after the expiration date.

10) Allergen panels are packed with desiccant and must be sealed after each use.

11) Some reagents included in the kit contain sodium azide as a preservative. Sodium azide has been reported to react with lead or copper plumbing to form potentially explosive metal azide, so always pour it along with a sufficient volume of water to prevent metal azide from forming in plumbing systems.

12) Substrate solution may contain black deposits. This is not abnormal and does not affect the test results.

13) Be careful to avoid the formation of bubbles. Especially when automatic devices are used, bubbles should be removed before starting the test, as bubbles can affect the distribution of solutions.

14) If the membranes are not completely dry after the final reaction, the test results analyzed by Q-SMART may be incorrect.

15) A product in a damaged package is not suitable for analysis.

[Packaging] Kit components, instructions for use.

[Storage instructions] Store at temperatures between 2°C and 8°C.

[Shelf life] Expiration date 18 months.

After opening, it is allowed to store at temperatures from 2 to 8 °C for no more than 3 months.

Transportation at temperatures from -15 °C to 30 °C for 5 days, by all types of transport in covered vehicles in accordance with the rules of transportation in force for this type of transport.

[Disposal]

When used in medical institutions, dispose of in accordance with class A and B SanPiN 2.1.7.3684-21. Dispose of this product in accordance with local, regional and national regulations.

[Manufacturer]

ProteomeTech Inc.



A-702, 401 Yangcheon-ro, Gangseo-gu, Seoul 07528, Korea.

[Authorized person]

On the territory of the Russian Federation, on issues related to the quality and circulation of a medical device, please contact:

109387, Moscow, Tikhaya str., 28, office 216; Limited Liability Company "BIO FOCUS" (LLC "BIO FOCUS"), biofocus@yandex.ru

[Symbols used on packaging and labels]

Symbol	Value	Symbol	Value
	Use up		Refer to instructions for use
	Batch code		In vitro diagnostic medical device
	Temperature Range from 2°C to 8°C		Manufacturer
			Authorised representative

	Caution! Refer to instructions for use		Number of panels / tests
	Not re-use		Biological risk
	Manufacturing date		

PROTIA *Allergy-Q* 64 Food panel

In vitro diagnostic kit for the quantitative determination of allergen-specific IgE antibodies in human serum or plasma using the immunoblotting method

[Appointment]

The reagent kit is intended for the quantitative determination of allergen-specific IgE antibodies in human serum or plasma using the immunoblotting method.

[Functional purpose]

Diagnostic aid. The use of the product has no population and demographic aspects.

[Application area]

Clinical laboratory diagnostics.

Medical device for in vitro diagnostic use only.

[Potential consumers of the product]

Medical facilities of any profile.

[Professional level of potential users]

Clinical laboratory diagnostics doctor, medical laboratory technician.

[Possible side effects]

Side effects have not been identified.

[Indications for use]

Allergy diagnostics.

[Contraindications for use]

Severe hemophilia.

The allergen panel contains components of animal origin:

- monoclonal antibodies against immunoglobulin E, Cat. No. 100105 (source of origin: mouse hybridomas, according to the documentation of the company Medix Biochemica (Finland));
- biotin-bovine serum albumin (source of origin: bovine serum albumin Cat. No. LY-0081 (bovine serum), according to the documentation of ProteomeTech Inc. (Republic of Korea)).

The sample diluent contains components of animal origin:

- gelatine Cat. No. G6650 (source of origin: bull flank - Type B, according to the documentation of SIGMA ALDRICH (USA)).

The antibody solution contains components of animal origin:

- biotin-conjugated mouse antibodies to human IgE, Cat. No. 100429/100106 (source of origin: mouse hybridomas, according to the documentation of Medix Biochemica (Finland));
- bovine serum albumin, Cat. No. LY-0081 (origin: bovine serum, according to the documentation of Boval Biosolutions LLC (USA)).

The enzyme solution contains components of animal origin:

- Bovine serum albumin, Cat. No. V5591 (source of origin: bovine serum), according to the documentation of the company Promega (USA));
- bovine serum albumin, Cat. No. LY-0081 (source of origin: bovine serum, according to the documentation of Boval Biosolutions LLC (USA)).

The components included in the kit are non-toxic and safe in the stated concentrations. Materials of animal origin are well fixed on the membranes of the test strip, placed in a plastic container with holes, which excludes contact with the human body when the medical device is used as intended.

Solutions, which include materials of animal origin, are well sealed with lids, when using solutions, a pipette and gloves are used, which exclude contact with the human body when using a medical device for its intended purpose.

Thus, there are no components in the medical device that are toxic, aggressive and dangerous when using the medical device for its intended purpose and subject to basic safety measures.

Sampling, processing, storage and handling of these materials in accordance with the documentation of the supplier / manufacturer of these materials, as well as general safety rules for working with biological materials.

[Description of the target analyte, information about its scientific validity, principle of action]**Operating principle**

PROTIA Allergy-Q is a multiplex diagnostic kit for allergy diagnostics based on the principle of enzyme immunoassay, which includes nitrocellulose membranes with different allergens adsorbed on regular interval lines, which allows you to test dozens of specific allergens in one test. PROTIA Allergy-Q uses a new parallel membrane technique and allows a wider variety of allergens to be assessed in a single test than other products using only one membrane. Allergen-specific IgE antibodies present in patient samples bind to antigens and are immobilized on the membrane even after the washing step.

Immobilized IgE antibodies bind to biotin antibodies to human IgE, and then biotin binds to alkaline phosphatase-conjugated lecnavidin. The color develops after the addition of the substrate in the last incubation with the enzyme, and the color intensity is analyzed with a color measuring device.

Description of the target analyte, scientific validity

Immunoblotting - a highly sensitive method for detecting proteins, based on a combination of electrophoresis and enzyme-linked immunosorbent assay (ELISA).

The pathogen antigens are separated by electrophoresis in polyacrylamide gel, then transferred from the gel to activated paper or nitrocellulose membrane and developed using ELISA. On these strips ("blots") the patient's serum is applied. Then, after incubation, the unbound antibodies of the patient are washed off and the serum against human immunoglobulins labeled with the enzyme is applied. The complex formed on the strip [antigen + patient's antibody + antibody against human Ig] is detected by adding a chromogenic substrate that changes color under the action of an enzyme.

Atopic allergy is a hypersensitive immunological condition mediated by immunoglobulin E (IgE) antibodies. Lymphocytic B cells, stimulated by a specific allergen, produce IgE antibodies to the antigen. IgE antibodies bind to receptors on the surface of mast cells or basophilic leukocytes through the Fc region. Subsequent binding of the allergen to cell-bound specific IgE causes cell degranulation and the release of vasoactive amines, resulting in smooth muscle contraction, itching, edema and trans-mucosal extracellular fluids. The most common clinical manifestations of this biological process are hay fever, asthma, dermatitis, urticaria, and anaphylactic shock.

[Kit Components]

Name of component	Characteristics	Quantity
1	2	
Allergen panel	Plastic cassette, with a plastic frame, containing a nitrocellulose membrane with adsorbed allergens (S-1, S-2, S-3 - standard lines; C - control line; from 1 to 60 - allergen lines)	20 panels
Sample diluent	A 20 ml polymer bottle with a blue screw cap containing 10 ml of blue-violet liquid	10ml x 1
Antibody solution	A 20 ml polymer bottle with a yellow screw cap, containing 10 ml of a yellowish or light yellow liquid	10ml x 1
Enzyme solution	A 20 ml polymer bottle with a red screw cap containing 10 ml of a colorless or yellowish liquid	10ml x 1
Substrate solution	A 20 ml polymer bottle with a black screw cap, containing 10 ml of yellow or yellowish liquid	10ml x 1
Washing solution 20x	A 20 ml polymer bottle with a white screw cap, containing 10 ml of concentrated colorless liquid	10ml x 1
Instructions for use		1

• Additional equipment required

1. Manual method: IMMUNOLOGICAL PHOTOMETRIC ANALYZER "Q-SMART" FOR IN VITRO DIAGNOSTICS ON "PROTIA" PANELS.

2. Semi-automatic method: "Q-PROCESSOR" AUTOMATED DEVICE FOR PROTIA PANEL PROCESSING FOR IN VITRO DIAGNOSTICS.

3. Automatic method: AUTOMATIC IMMUNOLOGICAL PHOTOMETRIC ANALYZER "Q-STATION ELITE" FOR IN VITRO DIAGNOSTICS ON PROTIA PANELS.

[Appearance]

1. Panel for allergens: a plastic panel in which a membrane is fixed in a plastic frame.
2. Sample diluent: blue-violet liquid.
3. Antibody solution: yellow or light yellow liquid.
4. Enzyme solution: colorless or light yellow liquid.
5. Substrate solution: light yellow liquid.
6. Washing solution 20x: colorless liquid.
7. Instructions for use.

[Analysis methods]

1. Preparation of reagents and samples

(1) Preparation of the cleaning solution.

Dilute the detergent solution 20x with deionized water 20 times before use.

(2) Sample preparation

The test uses human serum or plasma (with anticoagulants sodium citrate, heparin, EDTA). Before testing, remove blood cells or any solids by centrifugation for 10-15 minutes at 2500-3000 rpm. Hemolyzed or contaminated samples may lead to incorrect results; therefore, it is not allowed to use samples with a high lipid content and (or) with signs of hemolysis, and (or) with visible microbial growth. Store serum or plasma samples at $2 \sim 8^{\circ}\text{C}$ when used for a short period (within 2 weeks), and at -15°C or lower for longer storage. Thawed samples must be thoroughly mixed before testing. When working with thawed samples, it should be borne in mind that each thawing-freezing cycle leads to the destruction of detected antibodies, reducing the quality of the test material. In order to obtain reliable results of the study, it is recommended to use only a single freeze-thaw of the samples.

2. Analysis procedure

Manual method

1) Optimal conditions for the laboratory - conducting the experiment at $+25 \pm 2^{\circ}\text{C}$ and at $50 \pm 20\%$ humidity.

2) All reagents should be brought to room temperature 30 minutes prior to use and mix thoroughly immediately before use. Open the packaging of the allergen panel after adjusting it to room temperature.

3) Fully moisten the test membrane by adding 300 μL of diluted detergent solution to the allergen panel, shake the panel for 5 minutes, and then remove any remaining detergent solution (100 rpm recommended).

4) Pipette 250 μL of Sample Diluent into the Allergen Plate.

5) Pipette 50 μL of patient sample into the panel and incubate with shaking at room temperature for 45 minutes.

6) Remove residual solutions from the panel and flush the membrane twice with diluted detergent solution. With each wash, add 300 μL of diluted wash solution, incubate with shaking for 5 minutes and drain. The solution should not remain in the panel.

7) Dispense 250 μL of antibody solution to the panel and incubate with shaking at room temperature for 30 minutes.

8) Flush the membrane according to method (6).

9) Pipette 250 μL of enzyme solution into the panel and incubate with shaking until for 30 minutes.

10) Flush the membrane according to method (6).

11) Dispense 250 μL of substrate solution and incubate with shaking at room temperature in a dark room for 20 minutes.

12) Remove the substrate solution from the panel and rinse it with 250 µl deionized water.

13) Leave the panel on the table to dry naturally (make sure the membrane completely dry).

14) Insert the panels into the IMMUNOLOGICAL PHOTOMETRIC "Q-SMART" ANALYZER FOR IN VITRO DIAGNOSTICS ON THE "PROTIA" PANELS and evaluate the results. For more details, see the Manual for IMMUNOLOGICAL PHOTOMETRIC ANALYZER "Q-SMART" FOR DIAGNOSTICS IN VITRO ON THE PANELS "PROTIA".

* to incubate and shake the panel, use any orbital shaker with a rotation speed (oscillation frequency) of 100 rpm. There are no restrictions on the choice of the device model.

Semi-automatic method

1) AUTOMATED DEVICE "Q-PROCESSOR" FOR PROCESSING PANELS "PROTIA" FOR DIAGNOSTICS IN VITRO is equipped with panels and reagents for allergens.

2) Refer to the operation manual AUTOMATED "Q-PROCESSOR" DEVICE FOR PROTIA PANEL PROCESSING FOR IN VITRO DIAGNOSTICS and let the reactions progress.

3) Insert the panels into the Q-SMART analyzer and evaluate the results. For more details, see the Manual for the IMMUNOLOGICAL PHOTOMETRIC ANALYZER "Q-SMART" FOR DIAGNOSTICS IN VITRO ON THE PANELS "PROTIA".

Automatic method

1) ANALYZER AUTOMATIC IMMUNOLOGICAL PHOTOMETRIC "Q-STATION ELITE" FOR DIAGNOSTICS IN VITRO ON PANEL "PROTIA" is equipped with panels and reagents for allergens.

2) Refer to the operation manual. AUTOMATIC IMMUNOLOGICAL PHOTOMETRIC ANALYZER "Q-STATION ELITE" FOR DIAGNOSTICS IN VITRO ON THE PANELS "PROTIA" and give the reactions progress to eat.

3. Evaluation and interpretation of results

1) The amount of allergen-specific IgE antibodies is determined in IU / ml using the Q-SMART Analyzer or the Q-STATION ELITE Analyzer, the class is determined from the table below.

2) The amount of total IgE (tIgE) is expressed in IU / ml and can be quantified in the range from 1 to 2000 IU / ml.

Allergen-specific IgE		
ME / ml IgE	Class	Number
<0.15 ~ 0.34	0	Not found
0.35 ~ 0.69	1	Small
0.70 ~ 3.49	2	Moderate
3.50 ~ 17.49	3	Moderately large
17.50 ~ 49.99	4	Large
50.00 ~ 99.99	5	Very large
> 100	6	Extremely large

4. Quality control

The reference line (C) must be clearly visible. If the reference line is faintly visible, it cannot be read with optical measuring instruments. Retesting is recommended.

5. Characteristics

1) Analytical sensitivity (IgE detection threshold): 0.15 ME / ml.


2) Analytical specificity (interference and cross-reactions): None Interference and cross-reactivity of the kit with the following materials were found: bilirubin, albumin, immunoglobulins (IgA, IgM, IgG or IgD), 3 times higher than normal physiological levels.

3) Diagnostic sensitivity: 99.40% - 100% (with confidence probability 95%).

4) Diagnostic specificity: 98.77% - 100% (with a confidence level of 95%).

5) Reproducibility: 100%.

6) Coefficient of Variation (CV): less than 8%.

No	Allergen	The code	Group
1	Total IgE	tlgE	IgE
2	House dust	hi	General
3	<i>O. pteronyssinus</i>	d1	
4	<i>O. farinae</i>	d2	
6	Cat epithelium and dandruff	e1	
6	Dog dandruff	e5	
7	Egg white	f1	
8	Milk	f2	
9	Cockroach	i6	
10	Peanut	f13	
11	Soya beans	f14	General
12	Wheat flour	f4	
13	Alder	t2	
14	White birch	t3	
15	Oak	t7	
16	Ambrosia wormwood high	w1	
17	Hop	w22	
18	Common wormwood (Chernobyl)	w6	
19	<i>Alternaria alternata</i>	m6	
20	<i>Cladosporium herbarum</i>	m2	
21	<i>Aspergillus fumigatus</i>	m3	
22	Crabs	f23	
23	Shrimp	f24	
24	Mackerel	f206	
25	Cultured rye	g12	
26	Ccd	o214	
27	Peach	f95	
28	Apple	f49	
29	Sesame	f10	
30	Cod	f3	Seafood
31	Tuna / Salmon	f40 / f41	
32	Flounder / Anchovy (Hamsa) / Pollock	f254 / f313 /	

		f413		
33	Lobster (lobster) / Pacific squid	f80 / f258		
34	Acne	f264		
35	Mussels / Oysters / Clam / Scallop	f37 / f290 / f207 / f338		
36	Silkworm pupa	-		
37	Pork	f26		
38	Beef	f27		
39	Chicken meat	f83		
40	Mutton	f88		
41	Cheddar cheese	f81		
42	Barley flour	f6		
43	Rice	f9		
44	Buckwheat flour	f11		
45	Baker's yeast	f45		
46	Corn	f8		
47	Carrot	f31		
48	Potatoes	f35		
49	Garlic / Onion	f47 / f48		
50	Celery	f85		
51	Cucumber	f244		
52	Tomatoes	f25		
53	Orange	f33		
54	Strawberry	f44		
55	Kiwi / Mango / Banana	f84 / f91 / f92		
56	Chestnut	f299		
57	Walnut	f256		
58	Hazelnut	f17		
59	Almond / Pine nut / Sunflower	f20 / f253 / -		
60	Cocoa	f93		

Meat

Grains

Vegetables

Fruits

Nuts

Other

[
Pre-
cau-
tions
for
Users]
1)For in
vitro
diag-
nostic
use
only.
2)

PROTIA Allergy-Q can be used to aid in clinical diagnosis, the result is preliminary. The obtained analysis result should be correlated with the clinical picture. To clarify the diagnosis, a doctor's consultation is required.

- 3) False positive results can be obtained due to cross-reactivity of antigens tested with other antigens.
- 4) Do not smoke, eat or drink in areas where samples or tester are used new reagents.
- 5) Optimal conditions for a laboratory are to conduct an experiment at $+ 25 \pm 2$ °C and at $50 \pm 20\%$ humidity.
- 6) All specimens may include unknown infectious materials. When handling samples, wear disposable gloves and wash your hands after testing.
- 7) To avoid accidents, do not place needles, knives, or any other object that could cause harm when handling specimens and reagents.
- 8) All patient samples and components used in the kit must be considered as sources of biological hazard. They should be disposed of in accordance with the relevant guidelines.
- 9) Do not use kits after the expiration date.
- 10) Allergen panels are packed with desiccant and must be sealed after each use.
- 11) Some reagents included in the kit contain sodium azide as a preservative. Sodium azide has been re-

ported to react with lead or copper plumbing to form potentially explosive metal azide, so always pour it along with a sufficient volume of water to prevent metal azide from forming in plumbing systems.

12) Substrate solution may contain black deposits. This is not abnormal and does not affect the test results.

13) Be careful to avoid the formation of bubbles. Especially when automatic devices are used, bubbles should be removed before starting the test, as bubbles can affect the distribution of solutions.

14) If the membranes are not completely dry after the final reaction, the test results analyzed by Q-SMART may be incorrect.

15) A product in a damaged package is not suitable for analysis.

[Packaging] Kit components, instructions for use.

[Storage instructions] Store at temperatures between 2°C and 8°C.

[Shelf life] Expiration date 18 months.

After opening, it is allowed to store at temperatures from 2 to 8 °C for no more than 3 months.

Transportation at temperatures from -15 °C to 30 °C for 5 days, by all types of transport in covered vehicles in accordance with the rules of transportation in force for this type of transport.

[Disposal]

When used in medical institutions, dispose of in accordance with class A and B SanPiN 2.1.7.3684-21. Dispose of this product in accordance with local, regional and national regulations.

[Manufacturer]

ProteomeTech Inc. (A-702, 401 Yangcheon-ro, Gangseo-gu, Seoul 07528, Korea.

[Authorized person]

On the territory of the Russian Federation, on issues related to the quality and circulation of a medical device, please contact:

109387, Moscow, Tikhaya str., 28, office 216; Limited Liability Company "BIO FOCUS" (LLC "BIO FOCUS"), biofocus@yandex.ru

[Symbols used on packaging and labels]

Symbol	Value	Symbol	Value
	Use up		Refer to instructions for use
	Batch code		In vitro diagnostic medical device
	Temperature Range from 2°C to 8°C		Manufacturer
			Authorised representative
	Caution! Refer to instructions for use		Number of panels / tests

	No re-use		Biological risk
	Manufacturing date		

PROTIA *Allergy-Q* 64 Respiratory panel

In vitro diagnostic kit for the quantitative determination of allergen-specific IgE antibodies in human serum or plasma using the immunoblotting method

[Appointment]

The reagent kit is intended for the quantitative determination of allergen-specific IgE antibodies in human serum or plasma using the immunoblotting method.

[Functional purpose]

Diagnostic aid. The use of the product has no population and demographic aspects.

[Application area]

Clinical laboratory diagnostics.

Medical device for in vitro diagnostic use only.

[Potential consumers of the product]

Medical facilities of any profile.

[Professional level of potential users]

Clinical laboratory diagnostics doctor, medical laboratory technician.

[Possible side effects]

Side effects have not been identified.

[Indications for use]

Allergy diagnostics.

[Contraindications for use]

Severe hemophilia.

The allergen panel contains components of animal origin:

- monoclonal antibodies against immunoglobulin E, Cat. No. 100105 (source of origin: mouse hybridomas, according to the documentation of the company Medix Biochemica (Finland));
- biotin-bovine serum albumin (source of origin: bovine serum albumin Cat. No. LY-0081 (bovine serum), according to the documentation of ProteomeTech Inc. (Republic of Korea)).

The sample diluent contains components of animal origin:

- gelatine Cat. No. G6650 (source of origin: bull flank - Type B, according to the documentation of SIGMA ALDRICH (USA)).

The antibody solution contains components of animal origin:

- biotin-conjugated mouse antibodies to human IgE, Cat. No. 100429/100106 (source of origin: mouse hybridomas, according to the documentation of Medix Biochemica (Finland));

- bovine serum albumin, Cat. No. LY-0081 (origin: bovine serum, according to the documentation of Boval Biosolutions LLC (USA)).

The enzyme solution contains components of animal origin:

- Bovine serum albumin, Cat. No. V5591 (source of origin: bovine serum), according to the documentation of the company Promega (USA));

- bovine serum albumin, Cat. No. LY-0081 (source of origin: bovine serum, according to the documentation of Boval Biosolutions LLC (USA)).

The components included in the kit are non-toxic and safe in the stated concentrations. Materials of animal origin are well fixed on the membranes of the test strip, placed in a plastic container with holes, which excludes contact with the human body when the medical device is used as intended.

Solutions, which include materials of animal origin, are well sealed with lids, when using solutions, a pipette and gloves are used, which exclude contact with the human body when using a medical device for its intended purpose.

Thus, there are no components in the medical device that are toxic, aggressive and dangerous when using the medical device for its intended purpose and subject to basic safety measures.

Sampling, processing, storage and handling of these materials in accordance with the documentation of the supplier / manufacturer of these materials, as well as general safety rules for working with biological materials.

[Description of the target analyte, information about its scientific validity, principle of action]

Operating principle

PROTIA Allergy-Q is a multiplex diagnostic kit for allergy diagnostics based on the principle of enzyme immunoassay, which includes nitrocellulose membranes with different allergens adsorbed on regular interval lines, which allows you to test dozens of specific allergens in one test. PROTIA Allergy-Q uses a new parallel membrane technique and allows a wider variety of allergens to be assessed in a single test than other products using only one membrane. Allergen-specific IgE antibodies present in patient samples bind to antigens and are immobilized on the membrane even after the washing step.

Immobilized IgE antibodies bind to biotin antibodies to human IgE, and then biotin binds to alkaline phosphatase-conjugated lecnavidin. The color develops after the addition of the substrate in the last incubation with the enzyme, and the color intensity is analyzed with a color measuring device.

Description of the target analyte, scientific validity

Immunoblotting - a highly sensitive method for detecting proteins, based on a combination of electrophoresis and enzyme-linked immunosorbent assay (ELISA).

The pathogen antigens are separated by electrophoresis in polyacrylamide gel, then transferred from the gel to activated paper or nitrocellulose membrane and developed using ELISA. On these strips ("blots") the patient's serum is applied. Then, after incubation, the unbound antibodies of the patient are washed off and the serum against human immunoglobulins labeled with the enzyme is applied. The complex formed on the strip [antigen + patient's antibody + antibody against human Ig] is detected by adding a chromogenic substrate that changes color under the action of an enzyme.

Atopic allergy is a hypersensitive immunological condition mediated by immunoglobulin E (IgE) antibodies. Lymphocytic B cells, stimulated by a specific allergen, produce IgE antibodies to the antigen. IgE antibodies bind to receptors on the surface of mast cells or basophilic leukocytes through the Fc region. Subse-

quent binding of the allergen to cell-bound specific IgE causes cell degranulation and the release of vasoactive amines, resulting in smooth muscle contraction, itching, edema and trans-mucosal extracellular fluids. The most common clinical manifestations of this biological process are hay fever, asthma, dermatitis, urticaria, and anaphylactic shock.

[Kit Components]

Name of component	Characteristics	Quantity
1	2	
Allergen panel	Plastic cassette, with a plastic frame, containing a nitrocellulose membrane with adsorbed allergens (S-1, S-2, S-3 - standard lines; C - control line; from 1 to 60 - allergen lines)	20 panels
Sample diluent	A 20 ml polymer bottle with a blue screw cap containing 10 ml of blue-violet liquid	10ml x 1
Antibody solution	A 20 ml polymer bottle with a yellow screw cap, containing 10 ml of a yellowish or light yellow liquid	10ml x 1
Enzyme solution	A 20 ml polymer bottle with a red screw cap containing 10 ml of a colorless or yellowish liquid	10ml x 1
Substrate solution	A 20 ml polymer bottle with a black screw cap, containing 10 ml of yellow or yellowish liquid	10ml x 1
Washing solution 20x	A 20 ml polymer bottle with a white screw cap, containing 10 ml of concentrated colorless liquid	10ml x 1
Instructions for use		1

- Additional equipment required

1. Manual method: IMMUNOLOGICAL PHOTOMETRIC ANALYZER "Q-SMART" FOR IN VITRO DIAGNOSTICS ON "PROTIA" PANELS.
2. Semi-automatic method: "Q-PROCESSOR" AUTOMATED DEVICE FOR PROTIA PANEL PROCESSING FOR IN VITRO DIAGNOSTICS.
3. Automatic method: ANALYZER AUTOMATIC IMMUNOLOGICAL PHOTOMETRIC "Q-STATION ELITE" FOR DIAGNOSTICS IN VITRO ON PANEL "PROTIA".

[Appearance]

1. Panel for allergens: a plastic panel in which a membrane is fixed in a plastic frame.
2. Sample diluent: blue-violet liquid.
3. Antibody solution: yellow or light yellow liquid.
4. Enzyme solution: colorless or light yellow liquid.

5. Substrate solution: light yellow liquid.
6. Washing solution 20x: colorless liquid.
7. Instructions for use.

[Analysis methods]

1. Preparation of reagents and samples

(1) Preparation of the cleaning solution.

Dilute the detergent solution 20x with deionized water 20 times before use.

(2) Sample preparation

The test uses human serum or plasma (with anticoagulants sodium citrate, heparin, EDTA). Before testing, remove blood cells or any solids by centrifugation for 10-15 minutes at 2500-3000 rpm. Hemolyzed or contaminated samples may lead to incorrect results; therefore, it is not allowed to use samples with a high lipid content and (or) with signs of hemolysis, and (or) with visible microbial growth. Store serum or plasma samples at 2 ~ 8 °C when used for a short period (within 2 weeks), and at -15 °C or lower for longer storage. Thawed samples must be thoroughly mixed before testing. When working with thawed samples, it should be borne in mind that each thawing-freezing cycle leads to the destruction of detected antibodies, reducing the quality of the test material. In order to obtain reliable results of the study, it is recommended to use only a single freeze-thaw of the samples.

2. Analysis procedure

Manual method

1) Optimal conditions for the laboratory - conducting the experiment at $+ 25 \pm 2$ °C and at $50 \pm 20\%$ humidity.

2) All reagents should be brought to room temperature 30 minutes prior to use and mix thoroughly immediately before use. Open the packaging of the allergen panel after adjusting it to room temperature.

3) Fully moisten the test membrane by adding 300 µl of diluted detergent solution to the allergen panel, shake the panel for 5 minutes, and then remove any remaining detergent solution (100 rpm recommended).

4) Pipette 250 µL of Sample Diluent into the Allergen Plate.

5) Pipette 50 µL of patient sample into the panel and incubate with shaking at room temperature for 45 minutes.

6) Remove residual solutions from the panel and flush the membrane twice with diluted detergent solution. With each wash, add 300 µl of diluted wash solution, incubate with shaking for 5 minutes and drain. The solution should not remain in the panel.

7) Dispense 250 µl of antibody solution to the panel and incubate with shaking at room temperature for 30 minutes.

8) Flush the membrane according to method (6).

9) Pipette 250 µL of enzyme solution into the panel and incubate with shaking until for 30 minutes.

10) Flush the membrane according to method (6).

11) Dispense 250 µl of substrate solution and incubate with shaking at room temperature in a dark room for 20 minutes.

12) Remove the substrate solution from the panel and rinse it with 250 µl deionized water.

13) Leave the panel on the table to dry naturally (make sure the membrane completely dry).

14) Insert the panels into the IMMUNOLOGICAL PHOTOMETRIC "Q-SMART" ANALYZER FOR IN VITRO DIAGNOSTICS ON THE "PROTIA" PANELS and evaluate the results. For more details, see the Manual for IMMUNOLOGICAL PHOTOMETRIC ANALYZER "Q-SMART" FOR DIAGNOSTICS IN VITRO ON THE PANELS "PROTIA".

* to incubate and shake the panel, use any orbital shaker with a rotation speed (oscillation frequency) of

100 rpm. There are no restrictions on the choice of the device model.

Semi-automatic method

- 1) AUTOMATED DEVICE "Q-PROCESSOR" FOR PROCESSING PANELS "PROTIA" FOR DIAGNOSTICS IN VITRO is equipped with panels and reagents for allergens.
- 2) Refer to the operation manual AUTOMATED "Q-PROCESSOR" DEVICE FOR PROTIA PANEL PROCESSING FOR IN VITRO DIAGNOSTICS and let the reactions progress.
- 3) Insert the panels into the Q-SMART analyzer and evaluate the results. For more details, see the Manual for IMMUNOLOGICAL PHOTOMETRIC ANALYZER "Q-SMART" FOR DIAGNOSTICS IN VITRO ON THE PANELS "PROTIA".

Automatic method

- 1) ANALYZER AUTOMATIC IMMUNOLOGICAL PHOTOMETRIC "Q-STATION ELITE" FOR DIAGNOSTICS IN VITRO ON PANEL "PROTIA" is equipped with panels and reagents for allergens.
- 2) Refer to the operation manual. AUTOMATIC IMMUNOLOGICAL PHOTOMETRIC ANALYZER "Q-STATION ELITE" FOR DIAGNOSTICS IN VITRO ON THE PANELS "PROTIA" and give the reactions progress to eat.

3. Evaluation and interpretation of results

- 1) The amount of allergen-specific IgE antibodies is determined in IU / ml using the Q-SMART Analyzer or the Q-STATION ELITE Analyzer, the class is determined from the table below.
- 2) Total IgE (tIgE) is expressed in IU / ml and can be quantified in the range of 1 to 2000 IU / ml.


Allergen-specific IgE		
ME / ml IgE	Class	Number
<0.15 ~ 0.34	0	Not found
0.35 ~ 0.69	1	Small
0.70 ~ 3.49	2	Moderate
3.50 ~ 17.49	3	Moderately large
17.50 ~ 49.99	4	Large
50.00 ~ 99.99	5	Very large
> 100	6	Extremely large

4. Quality control

The reference line (C) must be clearly visible. If the reference line is faintly visible, it cannot be read with optical measuring instruments. Retesting is recommended.

5. Characteristics

- 1) Analytical sensitivity (IgE detection threshold): 0.15 ME / ml.
- 2) Analytical specificity (interference and cross-reactions): None Interference and cross-reactivity of the kit with the following materials were found: bilirubin, albumin, immunoglobulins (IgA, IgM, IgG or IgD), 3 times higher than normal physiological levels.
- 3) Diagnostic sensitivity: 99.40% - 100% (95% confidence level).
- 4) Diagnostic specificity: 98.77% - 100% (with a confidence level of 95%).
- 5) Reproducibility: 100%.
- 6) Coefficient of Variation (CV): less than 8%.

No.	Allergen	Code	Group
1	Total IgE	tlgE	IgE
2	House dust	h1	General
3	<i>O. pteronyssinus</i>	d1	
4	<i>D. farinae</i>	d2	
5	Cat epithelium and dandruff	e1	
6	Dog dandruff	e5	
7	Egg white	f1	
8	Milk	f2	
9	Cockroach	i6	
10	Peanut	f13	
11	Soya beans	f14	
12	Wheat	f4	
13	Alder	t2	
14	Birch	t3	
15	Oak	t7	
16	Ambrosia wormwood high	w1	
17	Hop	w22	
18	Common wormwood (Chernobyl)	w6	
19	<i>Alternaria alternata</i>	m6	In room
20	<i>Cladosporium herbarum</i>	m2	
21	<i>Aspergillus fumigatus</i>	m3	
22	Crabs	f23	
23	Shrimp	f24	
24	Mackerel	f206	
25	Cultured rye	g12	
26	Ccd	o214	
27	Peach	f95	
28	Apple	f49	
29	Sesame	f10	
30	<i>Acarus siro</i>	d70	In room
31	<i>Tyrophagus putrescentiae</i>	d72	
32	Fragrant spikelet / Hedgehog combined Common reed / Polevitsa	g1 / g3 / g7 / g9	Pollen
33	Fingered pig	g2	
34	Timofeevka	g6	
35	<i>Penicillium notatum</i>	m1	
36	<i>Candida albicans</i>	m5	Pollen
37	Common hazel	t4	
38	European olive	t9	
39	Sycamore	t1 1	
40	White willow	t12	

[	41	Poplar	t14	
	42	Ash white	t15	
	43	White pine	t16	
	44	Japanese cedar	t17	
	45	Acacia	t19	
	46	Nivyanik (popovnik)	w7	
	47	Dandelion	w8	
	48	Plantain	w9	
	49	Ash pan (potash pot, stone)	w11	
	50	Goldenrod	w12	
	51	Amaranth (thrown back)	w14	
	52	Latex	k82	
	53	Bee venom	i1	
	54	Wasp poison	i3	
	55	Mouse / Rat	e71 / e73	Animals
	56	Rabbit	e82	
	57	The guinea pig	e6	
	58	Sheep epithelium	e81	
	59	Hamster	e84	
	60	Horse	e3	

Precautions for Users]

- 1) For in vitro diagnostic use only.
- 2) PROTIA Alergy-Q can be used to aid in clinical diagnosis, the result is preliminary. The obtained analysis result should be correlated with the clinical picture. To clarify the diagnosis, a doctor's consultation is required.
- 3) False positive results can be obtained due to cross-reactivity of antigens tested with other antigens.
- 4) Do not smoke, eat or drink in areas where samples or tester are used new reagents.
- 5) Optimal conditions for a laboratory are to conduct an experiment at $+ 25 \pm 2$ °C and at $50 \pm 20\%$ humidity.
- 6) All specimens may include unknown infectious materials. When handling samples, wear disposable gloves and wash your hands after testing.
- 7) To avoid accidents, do not place needles, knives, or any other object that could cause harm when handling specimens and reagents.
- 8) All patient samples and components used in the kit must be considered as sources of biohazards. They should be disposed of in accordance with the relevant guidelines.
- 9) Do not use kits after the expiration date.
- 10) Allergen panels are packed with desiccant and must be sealed after each use.
- 11) Some reagents included in the kit contain sodium azide as a preservative. Sodium azide has been reported to react with lead or copper plumbing to form potentially explosive metal azide, so always pour it along with a sufficient volume of water to prevent metal azide from forming in plumbing systems.
- 12) Substrate solution may contain black deposits. This is not abnormal and does not affect the test results.
- 13) Be careful to avoid the formation of bubbles. Especially when automatic devices are used, bubbles should be removed before starting the test, as bubbles can affect the distribution of solutions.
- 14) If the membranes are not completely dry after the final reaction, the test results analyzed by Q-SMART may be incorrect.
- 15) A product in a damaged package is not suitable for analysis.

[Packaging] Kit components, instructions for use.

[Storage instructions] Store at temperatures between 2°C and 8°C.

[Shelf life] Expiration date 18 months.

After opening, it is allowed to store at temperatures from 2 to 8 °C for no more than 3 months.

Transportation at temperatures from -15 °C to 30 °C for 5 days, by all types of transport in covered vehicles in accordance with the rules of transportation in force for this type of transport.

[Disposal and disposal]

When used in medical institutions, dispose of in accordance with class A and B SanPiN 2.1.7.3684-21. Dispose of this product in accordance with local, regional and national regulations.

[Manufacturer]

ProteomeTech Inc.





A-702, 401 Yangcheon-ro, Gangseo-gu, Seoul 07528, Korea.

[Authorized person]

On the territory of the Russian Federation, on issues related to the quality and circulation of a medical device, please contact:

109387, Moscow, Tikhaya str., 28, office 216; Limited Liability Company "BIO FOCUS" (LLC "BIO FOCUS"), biofocus@yandex.ru

[Symbols used on packaging and labels]

Symbol	Value	Symbol	Value
	Use up		Refer to instructions for use
	Batch code		In vitro diagnostic medical device
	Temperature Range from 2°C to 8°C		Manufacturer
			Authorised representative
	Caution! Refer to instructions for use		Number of panels / tests
	No re-use		Biological risk
	Manufacturing date		

PROTIA *Allergy-Q* 96M panel

In vitro diagnostic kit for the quantitative determination of allergen-specific IgE antibodies in human serum or plasma using the immunoblotting method

[Appointment]

The reagent kit is intended for the quantitative determination of allergen-specific IgE antibodies in human serum or plasma using the immunoblotting method.

[Functional purpose]

Diagnostic aid. The use of the product has no population and demographic aspects.

[Application area]

Clinical laboratory diagnostics.

Medical device for in vitro diagnostic use only.

[Potential consumers of the product]

Medical facilities of any profile.

[Professional level of potential users]

Clinical laboratory diagnostics doctor, medical laboratory technician.

[Possible side effects]


Side effects have not been identified.

[Indications for use]

Allergy diagnostics.

[Contraindications for use]

Severe hemophilia.

The allergen panel contains components of animal origin:

- monoclonal antibodies against immunoglobulin E, Cat. No. 100105 (source of origin: mouse hybridomas, according to the documentation of the company Medix Biochemica (Finland));
- biotin-bovine serum albumin (source of origin: bovine serum albumin Cat. No. LY-0081 (bovine serum), according to the documentation of ProteomeTech Inc. (Republic of Korea)).

The sample diluent contains components of animal origin:

- gelatin Cat. No. G6650 (source of origin: bull flank - Type B, according to the documentation of SIGMA ALDRICH (USA)).

The antibody solution contains components of animal origin:

- biotin-conjugated mouse antibodies to human IgE, Cat. No. 100429/100106 (source of origin: mouse hybridomas, according to the documentation of Medix Biochemica (Finland));
- bovine serum albumin, Cat. No. LY-0081 (origin: bovine serum, according to the documentation of Boval Biosolutions LLC (USA)).

The enzyme solution contains components of animal origin:

- Bovine serum albumin, Cat. No. V5591 (source of origin: bovine serum), according to the documentation of the company Promega (USA));
- bovine serum albumin, Cat. No. LY-0081 (source of origin: bovine serum, according to the documentation of Boval Biosolutions LLC (USA)).

The components included in the kit are non-toxic and safe in the stated concentrations. Materials of animal origin are well fixed on the membranes of the test strip, placed in a plastic container with holes, which excludes contact with the human body when the medical device is used as intended.

Solutions, which include materials of animal origin, are well sealed with lids, when using solutions, a pipette and gloves are used, which exclude contact with the human body when using a medical device for its intended purpose.

Thus, there are no components in the medical device that are toxic, aggressive and dangerous when using the medical device for its intended purpose and subject to basic safety measures.

Sampling, processing, storage and handling of these materials in accordance with the documentation of the supplier / manufacturer of these materials, as well as general safety rules for working with biological materials.

[Description of the target analyte, information about its scientific validity, principle of action]

Operating principle

PROTIA Allergy-Q is a multiplex diagnostic kit for allergy diagnostics based on the principle of enzyme immunoassay, which includes nitrocellulose membranes with different allergens adsorbed on regular interval lines, which allows you to test dozens of specific allergens in one test. PROTIA Allergy-Q uses a new parallel membrane technique and allows a wider variety of allergens to be assessed in a single test than other products using only one membrane. Allergen-specific IgE antibodies present in patient samples bind to antigens and are immobilized on the membrane even after the washing step.

Immobilized IgE antibodies bind to biotin antibodies to human IgE, and then biotin binds to alkaline phosphatase-conjugated levanidin. The color develops after the addition of the substrate in the last incubation with the enzyme, and the color intensity is analyzed with a color measuring device.

Description of the target analyte, scientific validity

Immunoblotting - a highly sensitive method for detecting proteins, based on a combination of electrophoresis and enzyme-linked immunosorbent assay (ELISA).

The pathogen antigens are separated by electrophoresis in polyacrylamide gel, then transferred from the gel to activated paper or nitrocellulose membrane and developed using ELISA. On these strips ("blots") the patient's serum is applied. Then, after incubation, the unbound antibodies of the patient are washed off and the serum against human immunoglobulins labeled with the enzyme is applied. The complex formed on the strip [antigen + patient's antibody + antibody against human Ig] is detected by adding a chromogenic substrate that changes color under the action of an enzyme.

Atopic allergy is a hypersensitive immunological condition mediated by immunoglobulin E (IgE) antibodies. Lymphocytic B cells, stimulated by a specific allergen, produce IgE antibodies to the antigen. IgE antibodies bind to receptors on the surface of mast cells or basophilic leukocytes through the Fc region. Subsequent binding of the allergen to cell-bound specific IgE causes cell degranulation and the release of vasoactive amines, resulting in smooth muscle contraction, itching, edema and trans-mucosal extracellular fluids. The most common clinical manifestations of this biological process are hay fever, asthma, dermatitis, urticaria, and anaphylactic shock.

[Kit Components]

Name of component	Characteristics	Quantity
-------------------	-----------------	----------

1	2	
Allergen panel	Plastic cassette, with a plastic frame, containing a nitrocellulose membrane with adsorbed allergens (S-1, S-2, S-3 - standard lines; C1, C2 - control lines; from 1 to 91 - allergen lines)	20 panels
Sample diluent	A 20 ml polymer bottle with a blue screw cap containing 10 ml of blue-violet liquid	10ml x 1
Antibody solution	A 20 ml polymer bottle with a yellow screw cap, containing 10 ml of a yellowish or light yellow liquid	10ml x 1
Enzyme solution	A 20 ml polymer bottle with a red screw cap containing 10 ml of a colorless or yellowish liquid	10ml x 1
Substrate solution	A 20 ml polymer bottle with a black screw cap, containing 10 ml of yellow or yellowish liquid	10ml x 1
Washing solution 20x	A 20 ml polymer bottle with a white screw cap, containing 10 ml of concentrated colorless liquid	10ml x 1
Instructions for use		1

• Additional equipment required

1. Manual method: IMMUNOLOGICAL PHOTOMETRIC ANALYZER "Q-SMART" FOR IN VITRO DIAGNOSTICS ON "PROTIA" PANELS.

2. Semi-automatic method: "Q-PROCESSOR" AUTOMATED DEVICE FOR PROTIA PANEL PROCESSING FOR IN VITRO DIAGNOSTICS.

3. Automatic method: ANALYZER AUTOMATIC IMMUNOLOGICAL PHOTOMETRIC "Q-STATION ELITE" FOR IN VITRO DIAGNOSTICS ON PROTIA PANELS.

[Appearance]

1. Panel for allergens: a plastic panel in which a membrane is fixed in a plastikova frame.
2. Sample diluent: blue-violet liquid.
3. Antibody solution: yellow or light yellow liquid.
4. Enzyme solution: colorless or light yellow liquid.
5. Substrate solution: light yellow liquid.
6. Washing solution 20x: colorless liquid.
7. Instructions for use.

[Analysis methods]

1. Preparation of reagents and samples

(1) Preparation of the cleaning solution.

Dilute the detergent solution 20x with deionized water 20 times before use.

(2) Sample preparation

The test uses human serum or plasma (with anticoagulants sodium citrate, heparin, EDTA). Before testing,

remove blood cells or any solids by centrifugation for 10-15 minutes at 2500-3000 rpm. Hemolyzed or contaminated samples may lead to incorrect results; therefore, it is not allowed to use samples with a high lipid content and (or) with signs of hemolysis, and (or) with visible microbial growth. Store serum or plasma samples at 2 ~ 8 °C when used for a short period (within 2 weeks), and at -15 °C or lower for longer storage. Thawed samples must be thoroughly mixed before testing. When working with thawed samples, it should be borne in mind that each thawing-freezing cycle leads to the destruction of detected antibodies, reducing the quality of the test material. In order to obtain reliable results of the study, it is recommended to use only a single freeze-thaw of the samples.

2. Analysis procedure

Manual method

- 1) Optimal conditions for the laboratory - conducting the experiment at $+ 25 \pm 2$ °C and at $50 \pm 20\%$ humidity.
 - 2) All reagents should be brought to room temperature 30 minutes prior to use and mix thoroughly immediately before use. Open the packaging of the allergen panel after adjusting it to room temperature.
 - 3) Fully moisten the test membrane by adding 500 µl of diluted detergent solution to the allergen panel, shake the panel for 5 minutes, and then remove any remaining detergent solution (100 rpm recommended).
 - 4) Pipette 300 µL Sample Diluent into the Allergen Plate.
 - 5) Pipette 100 µL of patient sample into the panel and incubate with shaking at room temperature for 45 minutes.
 - 6) Remove residual solutions from the panel and flush the membrane twice with diluted detergent solution. With each wash, add 400 µl of diluted wash solution, incubate with shaking for 5 minutes and drain. The solution should not remain in the panel.
 - 7) Dispense 400 µl of antibody solution to the panel and incubate with shaking at room temperature for 30 minutes.
 - 8) Flush the membrane according to method (6).
 - 9) Pipette 400 µL of enzyme solution into the panel and incubate with shaking until for 30 minutes.
 - 10) Flush the membrane according to method (6).
 - 11) Dispense 400 µl of substrate solution and incubate with shaking at room temperature in a dark room for 20 minutes.
 - 12) Remove the substrate solution from the panel and rinse it with 400 µL deionized water.
 - 13) Leave the panel on the table to dry naturally (make sure membrane is completely dry).
 - 14) Insert the panels into the IMMUNOLOGICAL PHOTOMETRIC "Q-SMART" ANALYZER FOR IN VITRO DIAGNOSTICS ON THE "PROTIA" PANELS and evaluate the results. For more details, see the Manual for IMMUNOLOGICAL PHOTOMETRIC ANALYZER "Q-SMART" FOR DIAGNOSTICS IN VITRO ON THE PANELS "PROTIA".
- * to incubate and shake the panel, use any orbital shaker with a rotation speed (oscillation frequency) of 100 rpm. There are no restrictions on the choice of the device model.

Semi-automatic method

- 1) AUTOMATED DEVICE "Q-PROCESSOR" FOR PROCESSING PANELS "PROTIA" FOR DIAGNOSTICS IN VITRO is equipped with panels and reagents for allergens.
- 2) Refer to the operation manual AUTOMATED "Q-PROCESSOR" DEVICE FOR PROTIA PANEL PROCESSING FOR IN VITRO DIAGNOSTICS and let the reactions progress.
- 3) Insert the panels into the Q-SMART analyzer and evaluate the results. For more details, see the Manual for the IMMUNOLOGICAL PHOTOMETRIC ANALYZER "Q-SMART" FOR DIAGNOSTICS IN VITRO ON THE PANELS "PROTIA".

Automatic method

1) ANALYZER AUTOMATIC IMMUNOLOGICAL PHOTOMETRIC "Q-STATION ELITE" FOR DIAGNOSTICS IN VITRO ON PANEL "PROTIA" is equipped with panels and reagents for allergens.

2) Refer to the operation manual. AUTOMATIC IMMUNOLOGICAL PHOTOMETRIC ANALYZER "Q-STATION ELITE" FOR DIAGNOSTICS IN VITRO ON THE PANELS "PROTIA" and give the reactions progress to eat.

3. Evaluation and interpretation of results

1) The amount of allergen-specific IgE antibodies is determined in IU / ml using the Q-SMART Analyzer or the Q-STATION ELITE Analyzer, the class is determined from the table below.

2) Total IgE (tIgE) is expressed in IU / ml and can be quantified in the range of 1 to 2000 IU / ml.

Allergen-specific IgE		
ME / ml IgE	Class	Number
<0.15 ~ 0.34	0	Not found
0.35 ~ 0.69	1	Small
0.70 ~ 3.49	2	Moderate
3.50 ~ 17.49	3	Moderately large
17.50 ~ 49.99	4	Large
50.00 ~ 99.99	5	Very large
> 100	6	Extremely large

4. Quality control

The reference line (C) must be clearly visible. If the reference line is faintly visible, it cannot be read with optical measuring instruments. Retesting is recommended.

5. Characteristics

1) Analytical sensitivity (IgE detection threshold): 0.15 ME / ml.


2) Analytical specificity (interference and cross-reactions): None Interference and cross-reactivity of the kit with the following materials were found: bilirubin, albumin, immunoglobulins (IgA, IgM, IgG or IgD), 3 times higher than normal physiological levels.

3) Diagnostic sensitivity: 99.40% - 100% (95% confidence level).

4) Diagnostic specificity: 98.77% - 100% (with a confidence level of 95%).

5) Reproducibility: 100%.

6) Coefficient of Variation (CV): less than 8%.

No.	Allergen	The code	Group
1	Total IgE	tIgE	IgE
2	House dust	h1	General
3	<i>D. pteronyssinus</i>	d1	
4	<i>D. farinae</i>	d2	

5	Cat epithelium and dandruff	e1	
6	Dog dandruff	e5	
7	Egg white	f1	
8	Milk	f2	
9	Cockroach	i6	
10	Peanut	f13	
11	Soya beans	f14	
12	Wheat	f4	
13	Alder	t2	
14	Birch	t3	
15	Oak	t7	
16	Ambrosia high	w1	
17	Hop	w22	
18	Common wormwood (Chernobyl)	w6	
19	<i>Alternaria alternata</i>	m6	
20	<i>Cladosporium herbarum</i>	m2	
21	<i>Aspergillus fumigatus</i>	m3	
22	Crabs	f23	
23	Shrimp	f24	
24	Mackerel	f206	
25	Cultured rye	g12	
26	Ccd	o214	
27	Peach	f95	
28	Apple	f49	
29	Sesame	f10	
30	<i>Acarus siro</i>	d70	Indoors
31	<i>Tyrophagus putrescentiae</i>	d72	
32	Sweet spikelet / Hedgehog combined / Common reed / Polevitsa	g1 / g3 / g7 / g9	
33	Fingered pig	g2	
34	Timofeevka	g6	
35	<i>Penicillium notatum</i>	m1	
36	<i>Candida albicans</i>	m5	
37	Common hazel	t4	
38	European olive	t9	
39	Sycamore	t11	
40	White willow	t12	
41	Poplar	t14	Pollen
42	Ash white	t15	
43	White pine	t16	
44	Japanese cedar	t17	
45	Acacia	t19	
46	Nivyanik (popovnik)	w7	
47	Dandelion	w8	
48	Plantain	w9	
49	Ash pan (potash pot, stone)	w11	
50	Goldenrod	w12	
51	Amaranth (thrown back)	w14	

52	Latex	k82	Other
53	Bee venom	i1	Insects
54	Wasp poison	i3	
55	Mouse / Rat	e71 / e73	Animals
56	Rabbit	e82	
57	The guinea pig	e6	
58	Sheep epithelium	e81	
59	Hamster	e84	
60	Horse	e3	Seafood
61	Cod	f3	
62	Tuna / Salmon	f40 / f41	
63	Flounder / Anchovy (hamsa) / Saida	f254 / f313 / f413	
64	Lobster (lobster) / Pacific squid	f80 / f258	
65	Acne	f264	
66	Mussels / Oysters / Clam / Scallop	f37 / f290 / f207 / f338	Meat
67	Silkworm pupa	-	
68	Pork	f26	
69	Beef	f27	
70	Chicken meat	f83	
71	Mutton	f88	
72	Cheddar cheese	f81	Grains
73	Barley flour	f6	
74	Rice	f9	
75	Buckwheat flour	f11	Vegetables
76	Baker's yeast	f45	
77	Corn	f8	
78	Carrot	f31	Fruits
79	Potatoes	f35	
80	Garlic / Onion	f47 / f48	
81	Celery	f85	
82	Cucumber	f244	
83	Tomatoes	f25	
84	Orange	f33	Nuts
85	Strawberry	f44	
86	Kiwi / Mango / Banana	f84 / f91 /	
87	Chestnut	f299	Other
88	Walnut	f256	
89	Hazelnut	f17	
90	Almond / Pine nut / Sunflower	f20 / f253 /	Other
91	Cocoa	f93	

[Precautions for Users]

- 1) For in vitro diagnostic use only.
- 2) PROTIA Allergy-Q can be used to aid in clinical diagnosis, the result is preliminary. The obtained analysis result should be correlated with the clinical picture. To clarify the diagnosis, a doctor's consultation is required.

- 3) False positive results can be obtained due to cross-reactivity of antigens tested with other antigens.
- 4) Do not smoke, eat or drink in areas where samples or tester are used new reagents.
- 5) Optimal conditions for a laboratory are to conduct an experiment at $+ 25 \pm 2$ °C and at $50 \pm 20\%$ humidity.
- 6) All specimens may include unknown infectious materials. When handling samples, wear disposable gloves and wash your hands after testing.
- 7) To avoid accidents, do not place needles, knives, or any other object that could cause harm when handling specimens and reagents.
- 8) All patient samples and components used in the kit must be considered as sources of biohazards. They should be disposed of in accordance with the relevant guidelines.
- 9) Do not use kits after the expiration date.
- 10) Allergen panels are packed with desiccant and must be sealed after each use.
- 11) Some reagents included in the kit contain sodium azide as a preservative. Sodium azide has been reported to react with lead or copper plumbing to form potentially explosive metal azide, so always pour it along with a sufficient volume of water to prevent metal azide from forming in plumbing systems.
- 12) Substrate solution may contain black deposits. This is not abnormal and does not affect the test results.
- 13) Be careful to avoid the formation of bubbles. Especially when automatic devices are used, bubbles should be removed before starting the test, as bubbles can affect the distribution of solutions.
- 14) If the membranes are not completely dry after the final reaction, the test results analyzed by Q-SMART may be incorrect.
- 15) A product in a damaged package is not suitable for analysis.

[Packaging] Kit components, instructions for use.

[Storage instructions] Store at temperatures between 20°C and 80°C.

[Shelf life] Expiration date 18 months.

After opening, it is allowed to store at temperatures from 2 to 8 °C for no more than 3 months.

Transportation at temperatures from -15 °C to 30 °C for 5 days, by all types of transport in covered vehicles in accordance with the rules of transportation in force for this type of transport.

[Disposal]

When used in medical institutions, dispose of in accordance with class A and B SanPiN 2.1.7.3684-21. Dispose of this product in accordance with local, regional and national regulations.

[Manufacturer]

ProteomeTech Inc.

A-702, 401 Yangcheon-ro, Gangseo-gu, Seoul 07528, Korea.

[Authorized person]

On the territory of the Russian Federation, on issues related to the quality and circulation of a medical device, please contact:

109387, Moscow, Tikhaya str., 28, office 216; Limited Liability Company "BIO FOCUS" (LLC "BIO FOCUS"), biofocus@yandex.ru

[Symbols used on packaging and labels]

Symbol	Value	Symbol	Value
	Use up		Refer to instructions for use

	Batch code		In vitro diagnostic medical device
	Temperature Range from 2oC to 8oC		Manufacturer
			Authorised representative
	Caution! Refer to in- structions for use		Number of panels / tests
	No re-use		Biological risk
	Manufacturing date		

APPENDIX «Г»

Instructions for ensuring equipment grounding.

1. Connect the Device to a grounded outlet only.
2. The current carrying capacity of the earth circuit must be at least equivalent to that of the mains power supply circuit.
3. Connection should be made only to sockets with grounding contacts, which ensure the connection of the grounding circuit of the first and disconnection of the latter relative to the mains power supply circuit.
4. Any grounding failure inside or outside the device, as well as operating the device without grounding, can lead to dangerous consequences. Intentional disconnection of grounding is not allowed.
5. It is FORBIDDEN to use the Device with damaged outer insulation of the power cord.
6. Before starting operation, it is necessary to check the condition of the Device case. Cracks, mechanical damage to the case should not be observed.
7. Do not allow the power cord to be twisted or pinched in any way. Do not place the power cord where it can be damaged.
8. Do not use 3-pin to 2-pin adapters with a mains plug.

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Proteometech Inc.

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CEO: Kook Jin Lim

appeared before me and admitted said
principal's subscription to the attached

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Belong to Seoul Southern District Prosecutors' Office

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365 Sinwol-ro, Yangcheon-gu, Seoul, Korea
(Sinjeong-dong, chongsol Building 2F)

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This office has been authorized by the
Minister of Justice, the Republic of Korea
to act as Notary Public since Feb., 7, 2015
under Law No.9416

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