Medtronic

Medtronic Inc.
710 Medtronic Parkway
Minneapolis MN 55432-5604 USA
www.medtronic.com

[Signature] 3/17/2020

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Venous Inlet Port	3/8" (1.0 cm)
Access Port	Female Luer Port
Arterial Sample Port	Female Luer Port
Recirculation Port	1/4" (0.6 cm)
Gas Inlet Port	1/4" (0.6 cm)
Gas Outlet Port	3/8" (1.0 cm) non-barbed
Water Ports	1/2" (1.3 cm) quick disconnects
Size	Adult
Weight, kg	0.774
Dimensions (LxWxH), mm	100 x 159 x 217
Max gas flow rate	Initial Gas Transfer (O ₂) 403 ± 6 mL/min
Min gas flow rate	Initial Gas Transfer (CO ₂) 302 ± 5 mL/min
Heat exchanger efficiency ratio	0,46
Pressure range for the blood circuit	Blood side pressure drop 112 ± 4 mmHg Maximum: 750 mmHg (99,9 hPa)
Pressure range for the gas circuit	Gas pressure drop: ≤15mmHg

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
Possibility and methods of integration with other medical devices

Interacts with products for carrying out artificial blood circulation, presented on the scheme below.

Model	Name	Статус
TP	Temperature probe	Registered on the territory of Russia Federation (RZN from 2015/2448 29.12.2018)
61399401093	Affinity™ oxygenator with integrated CVR holder	Not a medical device. Outside the scope of supply

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Technical specifications of MD:

The tolerance is $\pm 5\%$ for all values, unless otherwise specified

Affinity™ NT membrane oxygenator with Balance biosurface and cardiotomy/venous reservoir, model 541B

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Specifications	Value
Membrane type	Microporous Polypropylene Hollow Fibers
Membrane Surface Area	2.5 m ²
Static Priming Volume	280 ml
Residual Blood Volume	149 ml
Minimum gas:blood ratio	0.5:1
Maximum gas:blood ratio	2:1
Recommended Blood Flow Rate	1-7 L/min
Maximum Water Side Pressure	206.8 kPa (30 psi)
Arterial Outlet Port	3/8" (1.0 cm)
Venous Inlet Port	3/8" (1.0 cm)
Access Port	Female Luer Port
Arterial Sample Port	Female Luer Port

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Recirculation Port	1/4" (0.6 cm)
Gas Inlet Port	1/4" (0.6 cm)
Gas Outlet Port	3/8" (1.0 cm) non-barbed
Water Ports	1/2" (1.3 cm) quick disconnects
Size	Adult
Weight, kg	0.774
Dimensions (LxWxH), mm	100 x 159 x 217
Max gas flow rate	Initial Gas Transfer (O ₂)
Min gas flow rate	395 ± 8 mL/min Initial Gas Transfer (CO ₂) 280 ± 8 mL/min
Heat exchanger efficiency ratio	0,46
Pressure range for the blood circuit	Blood side pressure drop 116 ± 5 mmHg Min Flow LTLblood = 7 mmHg Max Flow LTLblood = 95 mmHg
Pressure range for the gas circuit	Min Flow UTLgas = 3 mmHg Max Flow UTLgas = 23 mmHg

Cardiotomy/venous reservoir

Specifications	Value
Reservoir volume capacity	4000 mL ± 10 mL
Recommended blood flow rate	1-7 L/min
Maximum cardiotomy flow rate	6 L/min ± 0.1 L/min
Minimum operating level	200 mL ± 100 mL
Cardiotomy filtration	30 µm depth filter
Venous inlet screen	200 µm
Final reservoir screen	150 µm
Size	Adult
Positive/negative pressure relief valve	0 to 5 mmHg positive pressure in excess of 150 mmHg vacuum
Weight, kg	1.028
Dimensions (LxWxH), mm	240 x 193 x 413

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Affinity™ NT membrane oxygenator with Balance biosurface and cardiotomy/venous reservoir, model 541B

Item Description	Material Of Manufacture
Cap-MLL, White, NV	Virgin Polycarbonate or Polypropylene
Cap-MLL, Red, NV	Virgin Polycarbonate or Polypropylene
Assembly-Sleeve, Adaptor, Long	Virgin Polycarbonate Resin
Adhesive-Medithane, Fh-3911-A-Potting	Clear Urethane Elastomer
Adhesive-Medithane, Fh-3911-B- Potting	Clear Urethane Elastomer

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
USA Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

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[illegible]

## Affinity NT™

### Integrated Uncoated CVR/Oxygenator with Trillium™ Biosurface

#### 1. Model

541T

#### 2. Description

The Affinity NT™ integrated uncoated cardioplegic reservoir (CVR)/oxygenator with Trillium™ biosurface is a single-use gas exchange device with integral stainless steel heat exchanger and cardioplegic reservoir. The primary blood contact surfaces of the oxygenator are coated with Trillium biosurface (a polymer coating containing nonleaching heparin). Venous blood is collected and defoamed while cardioplegic blood is collected, defoamed, and filtered before mixing with the venous blood. The venous blood is pumped through the heat exchanger for temperature regulation and then outside the hollow fibers in the membrane compartment for gas exchange.

##### 2.1. Specifications

|                                             |                                                            |
|---------------------------------------------|------------------------------------------------------------|
| <b>Oxygenator</b>                           |                                                            |
| Membrane type                               | Microporous polypropylene hollow fibers                    |
| Membrane surface area                       | 2.5 m <sup>2</sup>                                         |
| Static priming volume                       | 270 mL                                                     |
| Residual blood volume <sup>a</sup>          | 143 mL                                                     |
| Minimum gas:blood ratio                     | 0.5:1                                                      |
| Maximum gas:blood ratio                     | 2:1                                                        |
| Recommended blood flow rate                 | 1 to 7 L/min                                               |
| Maximum water side pressure                 | 206.8 kPa (30 psi)                                         |
| Arterial outlet port                        | 1.0 cm (3/8 in)                                            |
| Venous inlet port                           | 1.0 cm (3/8 in)                                            |
| Access port                                 | Female luer port                                           |
| Arterial sample port                        | Female luer port                                           |
| Recirculation port                          | 0.6 cm (1/4 in)                                            |
| Gas inlet port                              | 0.6 cm (1/4 in)                                            |
| Gas outlet port                             | 1.0 cm (3/8 in) nonbarbed                                  |
| Water ports                                 | 1.3 cm (1/2 in) quick disconnects                          |
| Size                                        | Adult                                                      |
| Maximum transit temperature                 | 50°C (122°F)                                               |
| <sup>a</sup> Tested with water              |                                                            |
| <b>Cardioplegic reservoir</b>               |                                                            |
| Reservoir volume capacity <sup>a</sup>      | 4000 mL                                                    |
| Recommended blood flow rate                 | 1 to 7 L/min                                               |
| Maximum cardioplegic flow rate <sup>b</sup> | 6 L/min                                                    |
| Minimum operating level <sup>c</sup>        | 200 mL                                                     |
| Cardioplegic filtration                     | 30 µm depth filter                                         |
| Venous inlet screen                         | 200 µm                                                     |
| Final reservoir screen                      | 150 µm                                                     |
| Size                                        | Adult                                                      |
| Positive/negative pressure relief valve     | 0 to 5 mmHg positive pressure in excess of 150 mmHg vacuum |

<sup>a</sup> The tolerance of the gauge used for volume capacity measurement is ±10 mL.

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<sup>c</sup> The tolerance of the gauge used for minimum operating level measurement is ±100 mL.

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<sup>1</sup> Technology licensed under agreement from Stratus Therapeutics, Limited, United Kingdom.

## 2.2. Accessories available separately

| Model number | Description                                     |
|--------------|-------------------------------------------------|
| TP           | Temperature probe                               |
| 61309401003  | Affinity™ oxygenator with integrated CVR holder |

## 3. Indications for use

The Affinity HT integrated uncased CVR/oxygenator with Trilium bio-surface is intended to be used in an extracorporeal perfusion circuit to collect venous and cardiomy suctioned blood, cool or warm the blood, and oxygenate and remove carbon dioxide from the blood during routine cardiopulmonary bypass (CPB) procedures up to 6 hours in duration.

## 4. Contraindications

Use the device only as indicated.

## 5. Warnings

Read all warnings, precautions, and instructions for use carefully prior to use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.

- The fluid pathway is sterile and nonpyrogenic. Inspect each package and device prior to use. Do not use the device if the package is opened or damaged, the device is damaged, or the protective caps are not in place.
- If air bubbles and/or leaks are observed during priming and/or operation, these conditions may result in air embolism to the patient and/or fluid loss. The extracorporeal circuit must be continuously monitored. Do not use the device if these conditions are observed.
- Clear all gas emboli from the extracorporeal circuit before initiating bypass. Gas emboli are hazardous to the patient.
- In any extracorporeal circuit, the use of a prebypass and arterial line filter should be considered.
- Monitoring of preoxygenator circuit pressure should be considered.
- The blood phase pressure must be greater than the gas phase pressure at all times in the oxygenator.
- Do not exceed blood phase pressure of 750 mmHg.
- Do not exceed heat exchanger water phase pressure of 206.8 kPa (30 psi).
- Heat exchanger water temperature shall not exceed 42°C (107°F).
- Follow institution CPB protocol when using disinfectants in the heat exchanger during bypass. The integrity of the water path has been verified with hydrogen peroxide (330 ppm). Contact Medtronic for information regarding the use of additional disinfectants.
- Band all tubing connections for added protection against high fluid pressures.
- A replacement oxygenator should be readily available during perfusion.
- Do not obstruct gas outlet vents to prevent the gas side pressures from exceeding blood side pressures.
- The use of a blood level sensing mechanism during operation of this device is recommended.
- Ensure that the venous reservoir outlet is always positioned above the highest point in the membrane compartment of the Affinity HT oxygenator.
- Do not use a vacuum source in excess of 160 mmHg suction through this device.
- Do not occlude or obstruct the vent port of the CVR during operation.
- Only persons thoroughly trained in CPB procedures should use this device. Operation of each device requires constant supervision by qualified personnel for patient safety. Perfusion must be carefully and constantly monitored.
- Each device is intended for single use only. Do not reuse or resterilize. Sterilized using ethylene oxide.
- Attach tubing in such a manner as to prevent kinks or restrictions that may alter blood, water, or gas flow.
- Do not allow alcohol, alcohol-based fluids, anesthetic liquids (such as isoflurane), or corrosive solvents (such as acetone) to come into contact with the device as they may jeopardize the structural integrity.

## 6. Precautions

- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize any component of the system. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device which could result in patient injury, illness, or death.
- Use aseptic technique in all procedures.
- Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures. Adequate anticoagulation must be maintained before and during bypass.
- Additional warnings and precautions applicable to specific procedures can be found at appropriate places in the instructions for use.
- It is the responsibility of the user to dispose of the device in accordance with local regulations and hospital procedures.

- This product contains phthalates.

#### 7. Adverse effects

The following adverse effects are associated with the use of the product:

- Blood loss
- Coagulopathy
- Death
- Embolism
- Excessive blood component activation or thrombogenicity
- Extinction
- Hemolysis
- Hemolytic anemia
- Hypercarbia
- Hypotension
- Hypovolemia
- Hypoxia
- Infection
- Ischemia
- Neurological dysfunction
- Organ dysfunction

#### 8. Information for products with Trillium biosurface


The primary blood-contacting surfaces of the product are coated with Trillium biosurface. This coated surface enhances blood compatibility and provides a blood-contacting surface that is thromboresistant. Trillium biosurface contains nonleaching heparin derived from porcine intestinal mucosa<sup>2</sup>.

**Caution:** Product coated with Trillium biosurface is intended for single use only. Reesterilization may adversely affect the Trillium biosurface.

**Caution:** Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures.

#### 9. Instructions for use

##### 9.1. System setup

1. Cardiomyvenous reservoir
2. Oxygenator
3. Arterial filter
4. Prebypass filter

Figure 7

1. Remove the devices from the packaging carefully to ensure sterile fluid pathway.  
**Warning:** Ensure that aseptic technique is used during all stages of setup and use of this system.  
**Warning:** Before removing the devices from the packaging, inspect the packaging and product for damage. If the package or product is damaged, do not use as sterility of the device may have been compromised and/or performance may be affected.
2. Place the integrated CVR and oxygenator securely into its holder.  
**Note:** The Affinity NT oxygenator and CVR may be separated and utilized with separate holders. To separate the oxygenator from the CVR, turn the oxygenator approximately 1/8 of a turn in a counter-clockwise direction and place in their respective holders.
3. Remove the obturator cap from the positive/negative pressure relief valve.

<sup>2</sup> Nonleaching is defined as heparin at a level  $\leq 0.1$  IU/mL, as measured through clinically relevant infusion conditions.

4. Connect the 1.3 cm (1/2 in) water lines to the inlet and outlet ports on the heat exchanger. Start the heat exchanger water flow, and check for leaks from the water compartment to the blood compartment. Discard the oxygenator if water is present in blood compartment, and set up the standby oxygenator.
  5. Connect all blood and gas lines to their appropriate connections. Do not occlude tubing in pump raceway.
  6. Clamp the tubing at the oxygenator inlet, oxygenator arterial outlet, and oxygenator recirculation line outlet.
- Note: A CO<sub>2</sub> flush of the oxygenator is not required, but is recommended before priming.

### 9.2. Priming

Figure 8

1. Fill the CVR with 1.5 to 2 L of priming solution.
2. Remove all bubbles in the section of tubing from the CVR outlet to the clamp by the oxygenator inlet.
3. Remove the clamp from the oxygenator inlet, and allow solution to gravity prime the heat exchanger and fiber bundle of the oxygenator.

Note: Albumin has been reported to reduce the rare instance of transient high inlet pressure phenomenon in membrane oxygenators. Therefore, albumin in the prime solution should be considered.

### 9.3. Recirculation

Figure 9

1. After placing tubing into a properly occluded pump raceway, remove clamp from recirculation outlet per Figure 9 to allow recirculation between the oxygenator and venous reservoir.
2. Begin recirculation flow at 500 mL/min and gradually increase flow to 5 L/min.
3. Ensure that oxygenator is bubble-free.
4. Gradually reduce flow, stop the pump, and clamp the recirculation line.

#### 9.4. Priming A/V circuit

Figure 10

1. Add additional prime solution to fill the rest of the perfusion circuit if needed. Uncrimp the arterial outlet; slowly pump prime solution into the remaining circuit. Purge air from the arterial sample port and access port at this time.
2. The arterial filter should now be primed per manufacturer's instructions.
3. Inspect the entire circuit and all components for leaks.
4. Upon completion of priming and debubbling, gradually reduce blood flow and stop the pump, close all purge lines, and clamp the arterial, venous, and recirculation lines.
5. Ensure that the entire extracorporeal circuit is free of bubbles before initiating bypass.

#### 9.5. Initiation of bypass

**Warning:** The blood phase pressure must be greater than the gas phase pressure at all times.

**Note:** Priming solution may be prewarmed through the heat exchanger before initiation of bypass.

1. Check for adequate levels of anticoagulation before initiation of bypass.
2. Remove the arterial and venous clamps and gradually increase blood flow. Then initiate gas to blood flow rate of 1:1 with a minimum  $\text{FiO}_2$  of 80%.
3. Initiate water flow through the heat exchanger.

**Warning:** Do not allow air to enter the circuit.

#### 9.6. Operation during bypass

1. Arterial  $\text{PO}_2$  is controlled by varying the percent concentration of oxygen present in the ventilating gas.
  - To decrease the  $\text{PO}_2$ , decrease the amount of oxygen in the ventilating gas by decreasing the  $\text{FiO}_2$  on your oxygen blender.
  - To increase the  $\text{PO}_2$ , increase the amount of oxygen in the ventilating gas by increasing the  $\text{FiO}_2$  on your oxygen blender.
2.  $\text{PCO}_2$  is controlled primarily by varying total gas flow rate.
  - To decrease  $\text{PCO}_2$ , increase the total gas flow rate to increase the amount of  $\text{CO}_2$  removed.
  - To increase  $\text{PCO}_2$ , decrease the total gas flow rate to decrease the amount of  $\text{CO}_2$  removed.
3. Patient temperature is controlled by regulating the temperature of water flow into the heat exchanger.

**Note:** Care should be taken to minimize temperature gradient during cooling and rewarming phases of bypass.

#### 9.7. Termination of bypass

1. Terminate bypass as individual case and patient condition indicate.
2. Turn off gas flow before turning off blood flow.
3. Turn off water flow to heat exchanger.
4. After terminating bypass, continue to slowly recirculate through oxygenator using recirculation line and/or reconnected arterial venous loop.

#### 9.8. Emergency oxygenator replacement

A standby Affinity HT integrated uncoated CVR/oxygenator with Triflim bio-surface should always be available during perfusion.

1. Turn off water flow to the heat exchanger, clamp and remove water lines. Remove oxygen line from gas inlet port.
2. Clamp venous line at reservoir. Turn off blood pump and clamp arterial outlet.
3. Remove any monitoring/sample line from oxygenator.

4. Double clamp recirculation line, oxygenator inlet and outlet tubing, and out tubing between clamps, leaving adequate lengths for reconnection.
  5. To remove oxygenator from holder, first release the heat exchanger support base by pulling the release pin on the base, and slide the base down the mast.
  6. To remove the oxygenator from the CVR, turn the oxygenator approximately 1/8 of a turn in a counter-clockwise direction.
  7. Attach the replacement oxygenator to the CVR by inserting the positioning peg located on the gas cap into the coupling slots on the CVR, and turn clockwise until the oxygenator locks into position.
  8. Replace the heat exchanger support base onto the holder by sliding the base back up the mast until the heat exchanger rests securely in the horizontal channel of the base. Secure the support base by inserting the release pin of the support base into the holder mast.
  9. Reconnect the recirculation, oxygenator inlet, and outlet tubing. Band all connections.
  10. Connect water lines, unclamp, turn on water source, and inspect for water leaks.
  11. Connect oxygen line.
  12. Remove clamps from recirculation and oxygenator inlet tubing.
- Warning:** Do not unclamp arterial line or venous return line at this time.
13. With sufficient volume in the CVR, turn on blood pump and slowly fill the oxygenator.
  14. Increase blood flow to 5 L/min through recirculation line.
  15. Stop pump and clamp recirculation line.
  16. Ensure that the entire system is free of leaks and gas bubbles.
  17. Remove all clamps from venous and arterial lines and reinstate bypass.
  18. **Cardiotomy/venous reservoir vacuum-assisted venous drainage (VAVD)**  
(For use with the VAVD kit)


Figure 11.

1. Follow previously listed instructions for use for setup with the following modifications as listed here in steps 2 through 3.
2. Attach a calibrated means of monitoring negative pressure to either the venous line or the venous reservoir.
3. If a 1.0 cm (3/8 in) venous line is to be used, first attach a venous port reducing adapter to the 1.3 cm (1/2 in) venous inlet of the Affinity NT CVR reservoir. Connect the venous line to this reducing adapter.
4. Follow previously listed instructions for use for priming with the following additions as listed here in steps 5 through 10.
5. Replace the vented white luer caps on the manifold sampling line ports with nonvented luer caps. Place plugs ("dead-enders") on all unused 0.8 cm (1/4 in) and 1.0 cm (3/8 in) ports.
6. Prime, debubble, and complete preparation of the circuit as per standard procedure. Prime and zero the venous line pressure monitor.

7. Before instituting bypass, prepare the vacuum portion of the setup. Attach a vented vacuum line to the 0.6 cm (1/4 in) reservoir vent port and to the vacuum regulator with vapor trap.  
**Warning:** Clamp the quick prime line before application of vacuum.  
**Warning:** Do not allow the vapor trap to become completely filled during use. This may allow fluid to enter the vacuum controller, or may prevent the reservoir from being vented to atmosphere when no vacuum is applied.
8. To vent the vacuum line, ensure that the "Y" connector with side tubing is in a location that is easily accessible and leave open to atmosphere. This may be clamped and unclamped regularly throughout the procedure to apply or stop vacuum (Figure 11).  
**Warning:** Do not exceed vacuum levels greater than -80 mmHg as measured in the venous line.
9. Check to see if the vacuum seals are secure before initiating bypass.  
**Note:** Vacuum will result in increased flow through the sampling manifold line during use.
10. Use vacuum assistance as needed to institute and maintain cardiopulmonary bypass.  
**Warning:** Do not apply vacuum to the venous reservoir when there is no forward blood flow through the oxygenator. This applies to both arterial centrifugal and roller pumps (roller may not be occlusive in all positions). This will prevent air from being pulled across the membrane into the blood path by the reservoir vacuum.  
**Warning:** Do not allow the venous reservoir to become pressurized as this could obstruct venous drainage, force air retrograde into the patient, or cause air to enter the blood path of the oxygenator.

#### 11. Additional information

The following additional information and data are available upon request: sterilization method, materials of the blood pathway, plasma leakage across the semipermeable membrane, blood and gas pathway pressure drop at the maximum blood and gas flow rates, blood cell damage, particle release, air-handling capability, airflow characteristics, breakthrough volume, filtration efficiency, relevant tolerances for data in the Instructions for Use, protocol summaries, and test methods used to determine coverage and leaching.

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- A. This LIMITED WARRANTY provides the following assurances to the patient who receives an Affinity NT integrated uncoated CVR/oxygenator with Trillium bio-surface, hereafter referred to as "The Product".
  - (1) Should the Product fail to function within normal tolerances due to a defect in materials or workmanship prior to its Use-by date, Medtronic will, at its option: (a) issue a credit to the purchaser equal to the Purchase price, as defined in Subsection A(2), against the purchase of the replacement Product or (b) provide a functionally comparable replacement Product at no charge.
  - (2) As used herein, Purchase Price shall mean the lesser of the net invoiced price of the original, or current functionally comparable, or replacement Product.
- B. To qualify for the LIMITED WARRANTY, these conditions must be met:
  - (1) The Product must be used prior to its Use-by date.
  - (2) The unused portion of the Product must be returned to Medtronic and shall be the property of Medtronic.
  - (3) The Product must not have been altered or subjected to misuse, abuse, or accident.
  - (4) The Product must be used in accordance with the labeling and instructions for use provided with the Product.
- C. This LIMITED WARRANTY is limited to its express terms, in particular:
  - (1) Except as expressly provided by this LIMITED WARRANTY, Medtronic is not responsible for any direct, incidental, or consequential damages based on any defect, failure, or malfunction of the product, whether the claim is based on warranty, contract, tort, or otherwise.
  - (2) THIS LIMITED WARRANTY IS MADE ONLY TO THE PATIENT IN WHOM THE PRODUCT WAS USED. AS TO ALL OTHERS, MEDTRONIC MAKES NO WARRANTY, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE, WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOMER OR OTHERWISE. NO EXPRESS OR IMPLIED WARRANTY TO THE PATIENT SHALL EXTEND BEYOND THE PERIOD SPECIFIED IN A(1) ABOVE. THIS LIMITED WARRANTY SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.
  - (3) The exclusions and limitation set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid. This LIMITED WARRANTY gives the patient specific legal rights. The patient may also have other rights which vary from state to state.

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**Notes:** This LIMITED WARRANTY is provided by Medtronic, 710 Medtronic Parkway, Minneapolis, MN 55432. It applies only in the United States. Areas outside the United States should contact their local Medtronic representative for exact terms of the LIMITED WARRANTY.

**13. Important notice — limited warranty (for countries outside the US)**

A. This LIMITED WARRANTY provides assurance for the purchaser who receives an Affinity NT integrated uncoated CVR/oxysterator with Trillium biosurface, hereafter referred to as the "Product", that, should the Product fail to function to specification, Medtronic will issue a credit, equal to the original Product purchase price (but not to exceed the value of the replacement Product) against the purchase of any Medtronic replacement Product used for that patient.

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- (1) The Product must be used prior to its Use-by date.
- (2) The Product must be returned to Medtronic within 60 days after use and shall be the property of Medtronic.
- (3) The Product may not have been used for any other patient.

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- (1) In no event shall any replacement credit be granted where there is evidence of improper handling, improper implantation or material alteration of the replaced Product.
- (2) Medtronic is not responsible for any incidental or consequential damages based on any use, defect, or failure of the Product, whether the claim is based on warranty, contract, tort, or otherwise.

D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

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


Affinity NT oxygenator with Trillium biosurface

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


## Technical specifications of MD:

The tolerance is  $\pm 5\%$  for all values, unless otherwise specified

## Affinity NT oxygenator with Trillium biosurface

| Specifications                  | Value                                                         |
|---------------------------------|---------------------------------------------------------------|
| Membrane type                   | Microporous Polypropylene Hollow Fibers                       |
| Membrane Surface Area           | 2.5 m <sup>2</sup>                                            |
| Static Priming Volume           | 270 ml                                                        |
| Residual Blood Volume           | 143mL                                                         |
| Minimum gas:blood ratio         | 0.5:1                                                         |
| Maximum gas:blood ratio         | 2:1                                                           |
| Recommended Blood Flow Rate     | 1-7 L/min                                                     |
| Maximum Water Side Pressure     | 206.8 kPa (30 psi)                                            |
| Arterial Outlet Port            | 3/8" (1.0 cm)                                                 |
| Venous Inlet Port               | 3/8" (1.0 cm)                                                 |
| Access Port                     | Female Luer Port                                              |
| Arterial Sample Port            | Female Luer Port                                              |
| Recirculation Port              | 1/4" (0.6 cm)                                                 |
| Gas Inlet Port                  | 1/4" (0.6 cm)                                                 |
| Gas Outlet Port                 | 3/8" (1.0 cm) non-barbed                                      |
| Water Ports                     | 1/2" (1.3 cm) quick disconnects                               |
| Size                            | Adult                                                         |
| Weight, kg                      | 0.774                                                         |
| Dimensions (LxWxH), mm          | 100 x 159 x 217                                               |
| Max gas flow rate               | Initial Gas Transfer (O <sub>2</sub> )<br>403 $\pm$ 6 mL/min  |
| Min gas flow rate               | Initial Gas Transfer (CO <sub>2</sub> )<br>302 $\pm$ 5 mL/min |
| Heat exchanger efficiency ratio | 0,46                                                          |

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Pressure range for the blood circuit

Blood side pressure drop  
112  $\pm$  4 mmHg  
Maximum: 750 mmHg (99,9 hPa)

Pressure range for the gas circuit

Gas pressure drop:  
 $\leq 15$  mmHg

## FIBER

|                           |                           |
|---------------------------|---------------------------|
| Fiber Outer Diameter      | 289-299 microns           |
| Fiber Wall Thickness      | 29-33 microns             |
| Fiber Tensile Strength    | $\geq 175$ grams/filament |
| Fiber Elongation to Break | $\geq 50\%$               |

## Cardiotomy/venous reservoir

| Specifications                          | Value                                                      |
|-----------------------------------------|------------------------------------------------------------|
| Reservoir volume capacity               | 4000 mL $\pm$ 10 mL                                        |
| Recommended blood flow rate             | 1-7 L/min                                                  |
| Maximum cardiotomy flow rate            | 6 L/min $\pm$ 0.1 L/min                                    |
| Minimum operating level                 | 200 mL $\pm$ 100 mL                                        |
| Cardiotomy filtration                   | 30 $\mu$ m depth filter                                    |
| Venous inlet screen                     | 200 $\mu$ m                                                |
| Final reservoir screen                  | 150 $\mu$ m                                                |
| Size                                    | Adult                                                      |
| Positive/negative pressure relief valve | 0 to 5 mmHg positive pressure in excess of 150 mmHg vacuum |
| Weight, kg                              | 1.028                                                      |
| Dimensions (LxWxH), mm                  | 240 x 193 x 413                                            |

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|                                         |                                                                   |
|-----------------------------------------|-------------------------------------------------------------------|
| Cardiotomy/venous reservoir             |                                                                   |
| Component description                   | Composition                                                       |
| ADHESIVE-GEL,<br>181M, W/O FL. A        | Polyurethane oligomer/(Meth)<br>Acrylate<br>Monomer               |
| ASSEMBLY-SLEEVE,<br>ADAPTOR, LONG       | Virgin Polycarbonate Resin,                                       |
| SCREEN-VENOUS<br>TUBE, AFFINITY CVR     | Polyester                                                         |
| SCREEN-HARDSHELL,<br>AFFINITY CVR       | Frame - Cryolite<br>Screen - Polyester                            |
| PORT-VENOUS,<br>COVER, AFFINITY<br>CVR  | Clear Polycarbonate                                               |
| FOAM-POLY,<br>CARDIOUT,<br>AFFINITY CVR | medical grade polyester,<br>thermally reticulated<br>polyurethane |
| WEB-TUBULAR                             | Polypropylene                                                     |
| TUBE-SEMI RIGID,<br>CVR                 | Medical grade PVC                                                 |
| TUBERECIRCULATION,<br>AFFINITY CVR      | Medical grade PVC                                                 |
| ADHV-UV CURABLE                         | Polyurethane oligomer/(Meth)<br>Acrylate<br>Monomer               |
| FOAM-ANTI, A                            | Octamethylcyclotetrasiloxane<br>(Silicone)                        |
| FOAM-POLY,<br>VENOUS TUBE, NT           | medical grade Polyester<br>thermally reticulated<br>polyurethane  |
| FOAM-POLY, CARDIIN,<br>NT               | medical grade Polyester and<br>NT                                 |

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## Affinity NT™

### Oxygenator with Balance™ Biosurface and Cardiotomy/Venous Reservoir

#### 1. Models

|      |                                                                                |
|------|--------------------------------------------------------------------------------|
| 511B | Affinity NT oxygenator with Balance biosurface <sup>a</sup>                    |
| 541B | Affinity NT oxygenator with Balance biosurface and cardiotomy/venous reservoir |

<sup>a</sup> Technology licensed under agreement from Biotechnology, Limited, United Kingdom.

#### Accessories sold separately

|             |                                                 |
|-------------|-------------------------------------------------|
| TP          | Temperature probe                               |
| 61390401093 | Affinity™ oxygenator with integrated CVR holder |

#### 2. Description

##### Oxygenator

The Affinity NT oxygenator with Balance biosurface is a single-use, microporous, hollow-fiber, gas-exchange device with plasma-resistant fiber and integrated heat exchanger. The primary blood-contacting surfaces of the oxygenator are coated with Balance biosurface to reduce platelet activation and adhesion and preserve platelet function.

##### Cardiotomy/venous reservoir

The Affinity NT uncoated cardiotomy/venous reservoir (CVR) is a single-use device designed to collect and store blood during extracorporeal circulation. Venous blood is collected and deoxygenated while cardiotomy blood is collected, deoxygenated, and filtered before mixing with the venous blood.

##### 2.1. Specifications

|                                           |                                                            |
|-------------------------------------------|------------------------------------------------------------|
| <b>Oxygenator</b>                         |                                                            |
| Membrane type                             | Microporous polypropylene hollow fibers                    |
| Membrane surface area                     | 2.5 m <sup>2</sup>                                         |
| Static priming volume                     | 280 mL                                                     |
| Residual blood volume <sup>a</sup>        | 148 mL                                                     |
| Minimum gas:blood ratio                   | 0.5:1                                                      |
| Maximum gas:blood ratio                   | 2:1                                                        |
| Recommended blood flow rate               | 1 to 7 L/min                                               |
| Maximum water side pressure               | 206.8 kPa (30 psi)                                         |
| Arterial outlet port                      | 1.0 cm (3/8 in)                                            |
| Venous inlet port                         | 1.0 cm (3/8 in)                                            |
| Access port                               | Female Luer port                                           |
| Arterial sample port                      | Female Luer port                                           |
| Recirculation port                        | 0.6 cm (1/4 in)                                            |
| Gas inlet port                            | 0.6 cm (1/4 in)                                            |
| Gas outlet port                           | 1.0 cm (3/8 in) nonbarbed                                  |
| Water ports                               | 1.3 cm (1/2 in) quick disconnects                          |
| Size                                      | Adult                                                      |
| Maximum transit temperature               | 50°C (122°F)                                               |
| <sup>a</sup> Tested with saline           |                                                            |
| <b>Cardiotomy/venous reservoir</b>        |                                                            |
| Reservoir volume capacity <sup>a</sup>    | 4000 mL                                                    |
| Recommended blood flow rate               | 1 to 7 L/min                                               |
| Maximum cardiotomy flow rate <sup>b</sup> | 6 L/min                                                    |
| Minimum operating level <sup>c</sup>      | 200 mL                                                     |
| Cardiotomy filtration                     | 30 µm depth filter                                         |
| Venous inlet screen                       | 200 µm                                                     |
| Final reservoir screen                    | 150 µm                                                     |
| Size                                      | Adult                                                      |
| Positive/negative pressure relief valve   | 0 to 5 mmHg positive pressure in excess of 150 mmHg vacuum |

<sup>a</sup> The tolerance of the gauge used for volume capacity measurement is ±18 mL.

<sup>b</sup> The tolerance of the gauge used for flow rate measurement is ±0.1 L/min.

<sup>c</sup> The tolerance of the gauge used for minimum operating level measurement is ±108 mL.

Store the product at room temperature.

Upon written request, a list of all materials used in the construction of the oxygenator and CVR will be made available.

### 3. Indications for use

The Affinity NT oxygenator with Gelance bio surface is intended to be used in an extracorporeal perfusion circuit to oxygenate and remove carbon dioxide from the blood and to cool or warm the blood during routine cardiopulmonary bypass (CPB) procedures up to 6 hours in duration.

The Affinity NT cardiotomy/venous reservoir (CVR) is intended to be used in an extracorporeal perfusion circuit to collect venous and cardiotomy suctioned blood during routine CPB procedures up to 6 hours in duration.

### 4. Contraindications

Use the device only as indicated.

### 5. Warnings

Read all warnings, precautions, and instructions for use carefully before use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.

- The fluid pathway is sterile and nonpyrogenic. Inspect each package and device before use. Do not use the device if the package is opened or damaged, the device is damaged, or the protective caps are not in place.
- Continuously monitor the extracorporeal circuit for air bubbles or leaks. Do not use the device if air bubbles or leaks are observed during priming or operation as they may result in air embolism to the patient or fluid loss.
- Do not operate the oxygenator and CVR outside the recommended parameter ranges.
- Once an oxygenator is primed with blood, adequate anticoagulation should be maintained per institution CPB protocol and the blood pathway should be constantly recirculated at a low flow rate.
- After initiation of bypass, blood flow should circulate through the oxygenator and CVR at all times within the recommended blood flow range during CPB, except during emergency oxygenator replacement.
- Ensure that unused ports remain capped to prevent contamination and that the caps are tightened to avoid leakage.
- Clear all gas emboli from the extracorporeal circuit before initiating bypass. Gas emboli are hazardous to the patient.
- In any extracorporeal circuit, the use of a prebypass and arterial line filter is recommended.
- Monitoring of circuit pressure is recommended.
- The blood phase pressure must be greater than the gas phase pressure at all times in the oxygenator.
- Do not exceed blood pressure of 100 kPa (750 mmHg) in the blood pathway of the oxygenator.
- Do not exceed heat exchanger water phase pressure of 206.8 kPa (30 psi).
- Heat exchanger water temperature shall not exceed 42°C (107°F).
- Follow institution CPB protocol when using disinfectants in the heater/cooler during bypass. The integrity of the water path has been verified with hydrogen peroxide (330 ppm). Contact Medtronic for information regarding the use of additional disinfectants.
- Band all tubing connections for added protection against high fluid pressures.
- A replacement oxygenator should be readily available during CPB.
- Do not obstruct gas outlet vents to prevent the gas side pressures from exceeding blood side pressures.
- Gas transfer rates may change over time, and adjustments of FIO<sub>2</sub> and gas flow rates may be needed to achieve desired gas transfer performance.
- The use of a blood level sensing mechanism during operation of this device is recommended.
- Ensure that the CVR outlet is always positioned above the highest point in the membrane compartment of the oxygenator.
- Do not use a vacuum source in excess of 150 mmHg suction through this device.
- Do not occlude or obstruct the vent port of the CVR during operation (except when using VAVD).
- Only persons thoroughly trained in CPB procedures should use this device. Operation of each device requires constant supervision by qualified personnel for patient safety. CPB must be carefully and constantly monitored.
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize any component of the system. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Each device is sterilized using ethylene oxide.
- Attach tubing in such a manner as to prevent kinks or restrictions that may alter blood, water, or gas flow.
- Do not add alcohol, alcohol-based fluids, anesthetic liquids (such as isoflurane), or corrosive solvents (such as acetone) to come into contact with the device as they may jeopardize the structural integrity.

### 6. Precautions

- Use aseptic technique in all procedures.

Instructions for Use English 27

master:20180508, altDocMod:20180409, altDocMod2:20171221, altDocStyle\_graph:20180318, altDocStyle:20180308, templates\_global:20180215, master\_pageast:20180222, pageaspec:20180129, body\_extant:20171025, languageLootup:20171025, table:20180606, character\_bin:20170404, 3aB:20110211, 3aB\_sug:20110211

- Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures. Maintain adequate anticoagulation per institution CPB protocol.
- Additional warnings and precautions applicable to specific procedures can be found at appropriate places in the instructions for use.
- It is the responsibility of the user to handle and dispose of the device in accordance with applicable laws, regulations, and hospital procedures.
- Applicable to products with phthalates identified on the product label: This product contains phthalates.

#### 7. Adverse effects

The following known adverse effects are associated with the use of the product: blood loss, coagulopathy, death, embolism, excessive blood component activation or thrombogenicity, embarginization, hemolysis, hemolytic anemia, hypercarbia, hypotension, hypovolemia, hypoxia, infection, ischemia, neurological dysfunction, and organ dysfunction.

#### 8. Information for products with Balance biosurface

The primary blood-contacting surfaces of the product are coated with Balance biosurface to reduce platelet activation and adhesion and preserve platelet function. See Figure 8.

Caution: Product coated with Balance biosurface is intended for single use only. Reesterification may adversely affect the Balance biosurface.

Caution: Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures.

#### 9. Instructions for use

##### 9.1. System setup

1. Remove the devices from the packaging carefully to ensure sterile fluid pathway.

Warning: Use aseptic technique during all stages of setup and use of this system.

Warning: Before removing the devices from the packaging, inspect the packaging and product for damage. If the package or product is damaged, do not use as sterility of the device may have been compromised and/or performance may be affected.

2. Place the integrated CVR and oxygenator securely into its holder.

Note: The Affinity HT oxygenator and CVR may be separated and utilized with separate holders. To separate the oxygenator from the CVR, turn the oxygenator approximately 1/8 of a turn in a counter-clockwise direction and place in their respective holders.

3. Remove the saturator cap from the positive/negative pressure relief valve.

4. Connect the 1.3 cm (1/2 in) water lines to the inlet and outlet ports on the heat exchanger. Start the heat exchanger water flow, and check for leaks from the water compartment to the blood compartment. Discard the oxygenator if water is present in blood compartment, and set up the standby oxygenator.

5. Connect all blood and gas lines to their appropriate connections. Do not occlude tubing in pump raceway.

6. Clamp the tubing at the oxygenator inlet, oxygenator arterial outlet, and oxygenator recirculation line outlet.  
Note: A CO<sub>2</sub> flush of the oxygenator is not required, but is recommended before priming.

7. The sampling manifold should be attached so that the tubing segment with the one-way duckbill valve is positioned between the arterial sampling location and the manifold. Arterial blood should flow through the duckbill valve toward the manifold.

8. If using a recirculation line, connect the recirculation line to the recirculation port on the oxygenator. Connect the other end of the recirculation line to the recirculation port on the CVR.

9. Connect the quick-connect 1.3 cm (1/2 in) water lines to the inlet and outlet ports on the oxygenator. Start the water flow and check for leaks from the water compartment to the blood compartment before priming.  
Caution: Do not use the oxygenator if water is present in the blood compartment.

10. Ensure that the CVR is vented to atmosphere by removing the vent port cap. If using the CVR for VAVD (see Section 10), connect the vacuum line to the vent port.

##### 9.2. Priming and recirculation

Note: A CO<sub>2</sub> flush of the circuit before priming may improve the ease of priming.

1. Fill the CVR with the appropriate volume of priming solution per institution CPB protocol.
2. Remove all bubbles in the section of tubing (roller pump head/foot) or centrifugal pump (if used) from the CVR outlet to the oxygenator inlet.
3. Gravity prime the oxygenator or prime by forward pump flow per institution CPB protocol.
4. Gradually begin recirculation flow using a prebypass filter per institution CPB protocol.
5. Ensure that the oxygenator is free of bubbles including unused ports.

**Warning:** Tightly seal unused ports with the provided caps.

6. Add additional priming solution, if needed, to prime the rest of the extracorporeal circuit.
  7. Upon completion of priming and debubbling, gradually reduce the flow, stop the pump, close all purge lines, and clamp arterial, venous, and recirculation lines.
  8. Ensure that the entire extracorporeal circuit is free of bubbles before initiating bypass.
- Note:** Priming solution may be prewarmed using a heat exchanger before initiating bypass.

#### 9.3. Initiation of bypass

**Warning:** The blood phase pressure should be greater than the gas phase pressure at all times.

1. Ensure that the gas outlet port is not obstructed.
2. Check for adequate levels of anticoagulation before initiating bypass.
3. Remove the arterial and venous clamps, and gradually increase blood flow.
4. Initiate gas flow with appropriate gas settings per institution CPB protocol and clinical judgment specific to the patient and procedure.
5. Adjust water temperature to meet clinical demands.

#### 9.4. Operation during bypass

1. Arterial  $PO_2$  is controlled by varying the percent concentration of oxygen present in the ventilating gas.
    - a. To decrease the  $PO_2$ , decrease the percent of oxygen in the ventilating gas by adjusting the  $FI_{O_2}$  on the gas blender.
    - b. To increase the  $PO_2$ , increase the percent of oxygen in the ventilating gas by adjusting the  $FI_{O_2}$  on the gas blender.
  2.  $PCO_2$  is controlled primarily by varying the total gas flow rate.
    - a. To decrease  $PCO_2$ , increase the total gas flow rate to increase the amount of  $CO_2$  removed.
    - b. To increase  $PCO_2$ , decrease the total gas flow rate to decrease the amount of  $CO_2$  removed.
  3. Patient temperature is controlled by regulating the temperature of the water flow in the heat exchanger.
- Warning:** After initiation of bypass, blood flow should circulate through the oxygenator and CVR at all times within blood flow range during CPB, except during emergency oxygenator or CVR replacement.

#### 9.5. Termination of bypass

Terminate bypass per institution CPB protocol as indicated by individual case and patient condition.

#### 9.6. Emergency oxygenator replacement

A standby oxygenator should always be available during CPB.

1. Turn off water flow to the heat exchanger; clamp and remove the water lines. Remove the oxygen line from the gas inlet port.
  2. Clamp the venous line at the CVR. Turn off the blood pump and clamp the arterial outlet.
  3. Remove any monitoring/sample line from the oxygenator.
  4. Double clamp the recirculation line, oxygenator inlet and outlet tubing, and cut the tubing between clamps, leaving adequate lengths for reconnection.
  5. To remove the oxygenator from the holder, first release the heat exchanger support base by pulling the release pin on the base, and slide the base down the mast.
  6. To remove the oxygenator from the CVR, turn the oxygenator approximately 1/8 of a turn in a counter-clockwise direction.
  7. Attach the replacement oxygenator to the CVR by inserting the positioning pegs located on the gas cap into the coupling slots on the CVR, and turn clockwise until the oxygenator locks into position.
  8. Replace the heat exchanger support base onto the holder by sliding the base back up the mast until the heat exchanger rests securely in the horizontal channel of the base. Secure the support base by inserting the release pin of the support base into the holder mast.
  9. Reconnect the recirculation, oxygenator inlet, and outlet tubing. Bend all connections.
  10. Connect the water lines, unclamp, turn on the water source, and inspect for water leaks.
  11. Connect the oxygen line.
  12. Remove the clamps from the recirculation and oxygenator inlet tubing.
- Warning:** Do not unclamp the arterial line or the venous line at this time.
13. With sufficient volume in the CVR, turn on the blood pump and slowly fill the oxygenator.
  14. Increase blood flow through the recirculation line.
  15. Stop the pump and clamp the recirculation line.
  16. Ensure that the entire system is free of leaks and gas bubbles.

17. Remove all clamps from the venous and arterial lines and reinstate bypass.

#### 18. Cardiectomy/venous reservoir, vacuum-assisted venous drainage (VAVD)

##### 18.1. Warnings

- Do not use a vacuum regulator with maximum negative vacuum pressure greater than -20 kPa (-150 mmHg). Use of excessive negative pressure may lead to hemolysis.
- Do not occlude or obstruct the vent port of the CVR during operation (except when using VAVD).
- Draw only the intended drug amounts into the syringe that will be used for drug administration into the CVR during VAVD. Negative pressure in the CVR could draw more drug from the syringe into the CVR than intended and lead to excess drug delivery.
- All AVF shunt lines must be closed to atmosphere before stopping the pump or using low blood flow rates during VAVD. Failure to close the lines may allow air to be drawn into the blood side of the oxygenator from the filter.
- Appropriate occlusion of the arterial roller pump head is necessary when using VAVD. Failure to appropriately occlude the pump may allow air to be drawn into the blood side of the oxygenator from the filter.
- During the use of VAVD technique, use arterial filtration.
- Gradually return to atmospheric pressure when ending use of VAVD. A sudden change in pressure may result in turbulent blood flow inside the CVR.
- The line between the oxygenator and centrifugal pump (if used) must be clamped before stopping the pump during the use of VAVD. Failure to clamp off the arterial line may allow air to be drawn into the blood side of the oxygenator from the filter. Use of an arterial line one-way valve between the oxygenator and centrifugal pump is recommended.
- A controlled vacuum regulator is required during VAVD.
- A positive/negative pressure relief valve is required during VAVD (included, see Figure 3).

##### 18.2. Precautions

- Use of a pressure measurement device on the CVR and an additional negative pressure relief valve that operates at -13 kPa (-100 mmHg) should be considered.
- Use a vapor trap during the use of VAVD per institutional CPB protocol.
- Do not allow the vapor trap to become completely filled during use. This may allow fluid to enter the vacuum controller, or prevent the CVR from being vented to atmosphere when no vacuum is applied.
- Do not apply vacuum to the CVR when there is no forward blood flow through the oxygenator. This applies to arterial centrifugal and roller pumps (roller may not be occlusive in all positions). This will prevent air from being pulled across the membrane into the blood path by the CVR vacuum.

##### 18.3. VAVD instructions for use

1. Follow previously listed instructions for use for setup with the following modifications as listed here in steps 2 through 3.
2. Attach a calibrated means of monitoring negative pressure to either the venous line or the CVR.
3. Follow previously listed instructions for use for priming with the following additions as listed here in steps 5 through 9.
4. Replace the vented white luer caps on the manifold sampling line ports with nonvented luer caps. Place plugs ("dead-enders") on all unused 0.5 cm (1/4 in) and 1.0 cm (3/8 in) ports.
5. Prime, debubble, and complete preparation of the circuit as per standard procedure. Prime and zero the venous line pressure monitor.
6. Before instituting bypass, prepare the vacuum portion of the setup. Attach a vented vacuum line to the 0.5 cm (1/4 in) CVR vent port and to the vacuum regulator with vapor trap.  
Warning: Clamp the quick prime line before application of vacuum.  
Warning: Do not allow the vapor trap to become completely filled during use. This may allow fluid to enter the vacuum controller, or may prevent the CVR from being vented to atmosphere when no vacuum is applied.
7. To vent the vacuum line, ensure that the "Y" connector with side tubing is in a location that is easily accessible and leave open to atmosphere. This may be clamped and unclamped regularly throughout the procedure to apply or stop vacuum.  
Warning: Do not exceed maximum vacuum levels per institutional CPB protocol.
8. Check to see if the vacuum seals are secure before instituting bypass.  
Note: Vacuum will result in increased flow through the sampling manifold line during use.
9. Use vacuum assistance as needed to institute and maintain CPB.

**Warning:** Do not apply vacuum to the CVR when there is no forward blood flow through the oxygenator. This applies to both arterial centrifugal and roller pumps (roller may not be occlusive in all positions). This will prevent air from being pulled across the membrane into the blood path by the reservoir vacuum.

**Warning:** Do not allow the CVR to become pressurized as this could obstruct venous drainage, force air retrograde into the patient, or cause air to enter the blood path of the oxygenator.

#### 11. Additional information

The following additional information and data are available upon request: sterilization method, materials of the blood pathway, plasma leakage across the semipermeable membrane, blood and gas pathway pressure drop at the maximum blood and gas flow rates, blood cell damage, particle release, air-handling capability, antimicrobial characteristics, breakthrough volume, filtration efficiency, relevant tolerances for data in the Instructions for Use, protocol summaries, and test methods used to determine coverage and leaching.

#### 12. Important notice — limited warranty (US customers only)

A. This LIMITED WARRANTY provides the following assurances to the patient who receives an Affinity NT oxygenator with Balance biosurface or an Affinity NT oxygenator with Balance biosurface and cardiomy-venous reservoir, hereafter referred to as "The Product".

(1) Should the Product fail to function within normal tolerances due to a defect in materials or workmanship prior to its Use-by date, Medtronic will, at its option: (a) issue a credit to the purchaser equal to the Purchase price, as defined in Subsection A(2), against the purchase of the replacement Product or (b) provide a functionally comparable replacement Product at no charge.

(2) As used herein, Purchase Price shall mean the lesser of the net invoiced price of the original, or current functionally comparable, or replacement Product.

B. To qualify for the LIMITED WARRANTY, these conditions must be met:

(1) The Product must be used prior to its Use-by date.

(2) The unused portion of the Product must be returned to Medtronic and shall be the property of Medtronic.

(3) The Product must not have been altered or subjected to misuse, abuse, or accident.

(4) The Product must be used in accordance with the labeling and instructions for use provided with the Product.

C. This LIMITED WARRANTY is limited to its express terms. In particular:

(1) Except as expressly provided by this LIMITED WARRANTY, Medtronic is not responsible for any direct, incidental, or consequential damages based on any defect, failure, or malfunction of the product, whether the claim is based on warranty, contract, tort, or otherwise.

(2) THIS LIMITED WARRANTY IS MADE ONLY TO THE PATIENT (IN WHOM THE PRODUCT WAS USED, AS TO ALL OTHERS, MEDTRONIC MAKES NO WARRANTY, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE, WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOMER OR OTHERWISE. NO EXPRESS OR IMPLIED WARRANTY TO THE PATIENT SHALL EXTEND BEYOND THE PERIOD SPECIFIED IN A(1) ABOVE. THIS LIMITED WARRANTY SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.

(3) The exclusions and limitation set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid. This LIMITED WARRANTY gives the patient specific legal rights. The patient may also have other rights which vary from state to state.

(4) No person has any authority to bind Medtronic to any representation, condition, or warranty except this LIMITED WARRANTY.

**Note:** This LIMITED WARRANTY is provided by Medtronic, 710 Medtronic Parkway, Minneapolis, MN 55432. It applies only in the United States. Areas outside the United States should contact their local Medtronic representative for exact terms of the LIMITED WARRANTY.

#### 13. Important notice — limited warranty (for countries outside the US)

A. This LIMITED WARRANTY provides assurance for the purchaser who receives an Affinity NT oxygenator with Balance biosurface or an Affinity NT oxygenator with Balance biosurface and cardiomy-venous reservoir, hereafter referred to as the "Product", that, should the Product fail to function to specification, Medtronic will issue a credit, equal to the original Product purchase price (but not to exceed the value of the replacement Product) against the purchase of any Medtronic replacement Product used for that patient.

The Warnings contained in the product labeling are considered an integral part of this LIMITED WARRANTY. Contact your local Medtronic representative to obtain information on how to process a claim under this LIMITED WARRANTY.

- B. To qualify for the **LIMITED WARRANTY**, these conditions must be met:
- (1) The Product must be used prior to its Use-by date.
  - (2) The Product must be returned to Medtronic within 90 days after use and shall be the property of Medtronic.
  - (3) The Product may not have been used for any other patient.
- C. This **LIMITED WARRANTY** is limited to its express terms. In particular:
- (1) In no event shall any replacement credit be granted where there is evidence of improper handling, improper implantation or material alteration of the replaced Product.
  - (2) Medtronic is not responsible for any incidental or consequential damages based on any use, defect, or failure of the Product, whether the claim is based on warranty, contract, tort, or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this **LIMITED WARRANTY** is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid.

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
### Possibility and methods of integration with other medical devices

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Each Custom Circuit packaged for ease of use

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## Technical specifications of MD:

The tolerance is  $\pm 5\%$  for all values, unless otherwise specified

## Affinity NT oxygenator with Balance biosurface, including Affinity NT oxygenator with Balance biosurface

| Specifications                  | Value                                                                               |
|---------------------------------|-------------------------------------------------------------------------------------|
| Membrane type                   | Microporous Polypropylene Hollow Fibers                                             |
| Membrane Surface Area           | 2.5 m <sup>2</sup>                                                                  |
| Static Priming Volume           | 280 ml                                                                              |
| Residual Blood Volume           | 149 ml                                                                              |
| Minimum gas:blood ratio         | 0.5:1                                                                               |
| Maximum gas:blood ratio         | 2:1                                                                                 |
| Recommended Blood Flow Rate     | 1-7 L/min                                                                           |
| Maximum Water Side Pressure     | 206.8 kPa (30 psi)                                                                  |
| Arterial Outlet Port            | 3/8" (1.0 cm)                                                                       |
| Venous Inlet Port               | 3/8" (1.0 cm)                                                                       |
| Access Port                     | Female Luer Port                                                                    |
| Arterial Sample Port            | Female Luer Port                                                                    |
| Recirculation Port              | 1/4" (0.6 cm)                                                                       |
| Gas Inlet Port                  | 1/4" (0.6 cm)                                                                       |
| Gas Outlet Port                 | 3/8" (1.0 cm) non-barbed                                                            |
| Water Ports                     | 1/2" (1.3 cm) quick disconnects                                                     |
| Size                            | Adult                                                                               |
| Weight, kg                      | 0.774                                                                               |
| Dimensions (LxWxH), mm          | 100 x 159 x 217                                                                     |
| Max gas flow rate               | Initial Gas Transfer (O <sub>2</sub> )                                              |
| Min gas flow rate               | 395 $\pm$ 8 mL/min<br>Initial Gas Transfer (CO <sub>2</sub> )<br>280 $\pm$ 8 mL/min |
| Heat exchanger efficiency ratio | 0,46                                                                                |

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|                                   |                                                             |
|-----------------------------------|-------------------------------------------------------------|
| Cardiotomy/venous reservoir       |                                                             |
| Component description             | Composition                                                 |
| ADHESIVE-GEL, 181M, W/O FL. A     | Polyurethane oglegomer/(Meth) Acrylate Monomer              |
| ASSEMBLY-SLEEVE, ADAPTOR, LONG    | Virgin Polycarbonate Resin,                                 |
| SCREEN-VENOUS TUBE, AFFINITY CVR  | Polyester                                                   |
| SCREEN-HARDSHELL, AFFINITY CVR    | Frame - Cryolite<br>Screen - Polyester                      |
| PORT-VENOUS, COVER, AFFINITY CVR  | Clear Polycarbonate                                         |
| FOAM-POLY, CARDIOUT, AFFINITY CVR | medical grade polyester, thermally reticulated polyurethane |
| WEB-TUBULAR                       | Polypropylene                                               |
| TUBE-SEMI RIGID, CVR              | Medical grade PVC                                           |

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
Name: Harsh Dharamshi  
Title: Sr. Regulatory Affairs Specialist  
Signature: *[Signature]*  
Date: 9/26/2018

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**A Medtronic Custom-Designed Perfusion Circuit**  
*Partnering to make a better Custom Perfusion System for your patients*  
 Simple to complex • Pre-assembled circuit with components  
 of your choice • Delivered when you need it

**Each Custom Circuit  
 packaged for ease of use**

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



**Affinity NT oxygenator with Trillium biosurface**

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


Cardiotomy/venous reservoir

**Schemes of medical device:**  
**Affinity NT oxygenator with Trillium biosurface**

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

- |                        |                                 |
|------------------------|---------------------------------|
| 1 – Gas Inlet          | 7 – Arterial Sample Port        |
| 2 – Gas Outlet         | 8 – Arterial Blood Outlet       |
| 3 – Water In           | 9 – Temperature Monitoring Port |
| 4 – Water Out          | 10 – Gas Vent                   |
| 5 – Access Port        | 11 – Heat Exchanger             |
| 6 – Recirculation Port | 12 – Venous Blood Inlet         |


**Affinity NT oxygenator with Balance biosurface, including:**

**Affinity NT oxygenator with Balance biosurface**

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

- |                        |                                 |
|------------------------|---------------------------------|
| 1 – Gas Inlet          | 7 – Arterial Sample Port        |
| 2 – Gas Outlet         | 8 – Arterial Blood Outlet       |
| 3 – Water In           | 9 – Temperature Monitoring Port |
| 4 – Water Out          | 10 – Gas Vent                   |
| 5 – Access Port        | 11 – Heat Exchanger             |
| 6 – Recirculation Port | 12 – Venous Blood Inlet         |

Cardiotomy/venous reservoir  
Side view

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

1 – Venous return inlet - 1.3 cm (1/2 in), rotatable

2 – Venous temperature monitoring adapter

3 – Recirculation port with tubing - 0.6 cm (1/4 in)

4 – Vent port - 0.6 cm (1/4 in), nonbarbed

5 – Luer lock, unfiltered

6 – Auxiliary cardiotomy reservoir inlet - 1.0 cm (3/8 in), unfiltered

7 – Quick prime - 0.6 cm (1/4 in), unfiltered

8 – Cardiotomy turret, rotatable 180°

9 – Luer lock, filtered - 2 each

10 – Quick prime - 0.6 cm (1/4 in), filtered

11 – Cardiotomy inlets - 4 x 0.6 cm (1/4 in), 1 x 1.0 cm (3/8 in)

17 – Arterial temperature monitoring adaptor

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|                                      |                                                                                     |
|--------------------------------------|-------------------------------------------------------------------------------------|
| Membrane type                        | Microporous Polypropylene Hollow Fibers                                             |
| Membrane Surface Area                | 2.5 m <sup>2</sup>                                                                  |
| Static Priming Volume                | 270 ml                                                                              |
| Residual Blood Volume                | 143mL                                                                               |
| Minimum gas:blood ratio              | 0.5:1                                                                               |
| Maximum gas:blood ratio              | 2:1                                                                                 |
| Recommended Blood Flow Rate          | 1-7 L/min                                                                           |
| Maximum Water Side Pressure          | 206.8 kPa (30 psi)                                                                  |
| Arterial Outlet Port                 | 3/8" (1.0 cm)                                                                       |
| Venous Inlet Port                    | 3/8" (1.0 cm)                                                                       |
| Access Port                          | Female Luer Port                                                                    |
| Arterial Sample Port                 | Female Luer Port                                                                    |
| Recirculation Port                   | 1/4" (0.6 cm)                                                                       |
| Gas Inlet Port                       | 1/4" (0.6 cm)                                                                       |
| Gas Outlet Port                      | 3/8" (1.0 cm) non-barbed                                                            |
| Water Ports                          | 1/2" (1.3 cm) quick disconnects                                                     |
| Size                                 | Adult                                                                               |
| Weight, kg                           | 0.774                                                                               |
| Dimensions (LxWxH), mm               | 100 x 159 x 217                                                                     |
| Max gas flow rate                    | Initial Gas Transfer (O <sub>2</sub> )                                              |
| Min gas flow rate                    | 403 $\pm$ 6 mL/min<br>Initial Gas Transfer (CO <sub>2</sub> )<br>302 $\pm$ 5 mL/min |
| Heat exchanger efficiency ratio      | 0,46                                                                                |
| Pressure range for the blood circuit | Blood side pressure drop<br>112 $\pm$ 4 mmHg<br>Maximum: 750 mmHg (99,9 hPa)        |
| Pressure range for the gas circuit   | Gas pressure drop:<br>$\leq 15$ mmHg                                                |

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#### Affinity NT oxygenator with Balance biosurface, including

Affinity NT oxygenator with Balance biosurface

| Specifications              | Value                                   |
|-----------------------------|-----------------------------------------|
| Membrane type               | Microporous Polypropylene Hollow Fibers |
| Membrane Surface Area       | 2.5 m <sup>2</sup>                      |
| Static Priming Volume       | 280 ml                                  |
| Residual Blood Volume       | 149 ml                                  |
| Minimum gas:blood ratio     | 0.5:1                                   |
| Maximum gas:blood ratio     | 2:1                                     |
| Recommended Blood Flow Rate | 1-7 L/min                               |
| Maximum Water Side Pressure | 206.8 kPa (30 psi)                      |
| Arterial Outlet Port        | 3/8" (1.0 cm)                           |
| Venous Inlet Port           | 3/8" (1.0 cm)                           |
| Access Port                 | Female Luer Port                        |
| Arterial Sample Port        | Female Luer Port                        |
| Recirculation Port          | 1/4" (0.6 cm)                           |
| Gas Inlet Port              | 1/4" (0.6 cm)                           |

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|                                      |                                                                                                       |
|--------------------------------------|-------------------------------------------------------------------------------------------------------|
| Gas Outlet Port                      | 3/8" (1.0 cm) non-barbed                                                                              |
| Water Ports                          | 1/2" (1.3 cm) quick disconnects                                                                       |
| Size                                 | Adult                                                                                                 |
| Weight, kg                           | 0.774                                                                                                 |
| Dimensions (LxWxH), mm               | 100 x 159 x 217                                                                                       |
| Max gas flow rate                    | Initial Gas Transfer (O <sub>2</sub> )                                                                |
| Min gas flow rate                    | 395 ± 8 mL/min<br>Initial Gas Transfer (CO <sub>2</sub> )<br>280 ± 8 mL/min                           |
| Heat exchanger efficiency ratio      | 0,46                                                                                                  |
| Pressure range for the blood circuit | Blood side pressure drop<br>116 ± 5 mmHg<br>Min Flow LTLblood = 7 mmHg<br>Max Flow LTLblood = 95 mmHg |
| Pressure range for the gas circuit   | Min Flow UTLgas = 3 mmHg<br>Max Flow UTLgas = 23 mmHg                                                 |

#### Cardiotomy/venous reservoir

| Specifications                          | Value                                                      |
|-----------------------------------------|------------------------------------------------------------|
| Reservoir volume capacity               | 4000 mL ± 10 mL                                            |
| Recommended blood flow rate             | 1-7 L/min                                                  |
| Maximum cardiotomy flow rate            | 6 L/min ± 0.1 L/min                                        |
| Minimum operating level                 | 200 mL ± 100 mL                                            |
| Cardiotomy filtration                   | 30 µm depth filter                                         |
| Venous inlet screen                     | 200 µm                                                     |
| Final reservoir screen                  | 150 µm                                                     |
| Size                                    | Adult                                                      |
| Positive/negative pressure relief valve | 0 to 5 mmHg positive pressure in excess of 150 mmHg vacuum |
| Weight, kg                              | 1.028                                                      |
| Dimensions (LxWxH), mm                  | 240 x 193 x 413                                            |

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|                                                       |                                                                                                                                                 |
|-------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Affinity NT oxygenator with Balance biosurface</b> |                                                                                                                                                 |
| Item Description                                      | Material Of Manufacture                                                                                                                         |
| Cap-MLL, White, NV                                    | 100% Virgin Polycarbonate Type 104 or Polypropylene, Hercules Type PD626, or Flint Hills Resources P4G4Z-011 With Colorant, Ferro Type CH1010SP |
| Cap-MLL, Red, NV                                      | 100% Virgin Polycarbonate Type 104 or Polypropylene, Hercules Type PD626, or Flint Hills Resources P4G4Z-011 With Colorant, Meyer 2-1785        |
| Assembly-Sleeve, Adaptor, Long                        | 100% Virgin Polycarbonate Resin, Makrolon Rx1452-550115 Clear                                                                                   |
| Adhesive-Medithane, Fh-3911-A-Potting                 | (Clear Urethane Elastomer) Medithane FH-3911-A                                                                                                  |
| Adhesive-Medithane, Fh-3911-B- Potting                | (Clear Urethane Elastomer) Medithane FH-3911-B                                                                                                  |

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\*Technology licensed under agreement from BioInteractions, United, United Kingdom.

Barcode Area  
FOR PLACEMENT ONLY

Medtronic, Inc.  
710 Medtronic Parkway  
Minneapolis, MN 55432 USA

CE  
0123

M50605A001 Rev 1A

Example of Unit Label

# Medtronic Affinity NT™

REF 541T  
Catalog number

an Integrated Uncoated CVR/Oxygenator with Trillium™ Biosurface\* /  
bg интегриран CVR без покритие/Оксигенатор с биопокрытие Trillium™ /  
es Nepotlažený integrovaný kardiotorakální/venózní rezervoár (CVR)/Oxygenátor  
s hliníkovým Trillium™ / de Integrierter CVR ohne Beschichtung/Oxygenator mit  
Trillium™ bioverfläche / de Oxygenator mit integriertem unbeschichtetem CVR mit  
Trillium™ Biosurface / el Οξυγονωτής με βιοσφαιρικό Trillium™, με  
ενσωματωμένη μη επικαλυπτική δεξιά CVR / es Oxigenador con superficie  
biocompatible Trillium™ con reservorio de cardiotoraxial/venoso integrado no  
revestido / et Integreeritud katmata CVR / Trillium™ Biosurface kaaliga  
oksigenaator / fi Integroitu pinnottomaton kardiotoraxia- / venosivarsäiliö ja  
Trillium™-biopinnotteella päällystetty hapetin / fr Oxigénéateur avec réservoir de  
cardiotoraxio-veneux non enduit intégré et biosurface Trillium™ / hr Integrirani,  
neobloženi CVR/oksigenator s biopovršinom Trillium™ / hu Beépített, bevonat  
nélküli kardiotoraxiás vagy vénás tartály / oxigenátor Trillium™ biofelülettel /  
it Ossigenatore con cardiotoraco/venoso integrato non rivestito, con  
superficie biopassiva Trillium™ / it Integruotas nedengtas CVR / oksigenatorius su  
„Trillium™“ biologiškai dengtas / iv Integreliate kardiotorakali/venozais rezervuārs bez  
pārklājuma/oksigenators ar Trillium™ bioloģisko pārklājumu / nl Oxygenator met  
geïntegreerd ongecoat CVR en Trillium™-bio-oppervlakte/aag / no Integrert  
ubelagt CVR/oksigenator med Trillium™ biooverflate / pt Oxigenador  
z biopowłoką Trillium™ zintegrowany z niepowlekanym zbiornikiem kardiotorakalnym/  
złoty (CVR) / pt-br Oxigenador com biossuperfície Trillium™ e CVR integrado  
não revestido / ro Oxigenator cu biosuprafață Trillium™ și rezervor de  
cardiotoraxio-venos (CVR) fără înveliș, integrat / ru Интегрированный КБР без  
покрытия/оксигенатор с биопокрытием Trillium™ / sk Oxigenátor  
s integrovaným kardiotorakálnym/venóznym zásuvníkom bez povlaku  
a s biologickou povrchovou úpravou Trillium™ Biosurface / el Oksigenator z  
integriranim neprevlačenim rezervoarjem za venosko in kardiotorakalno kri z  
biološko površino Trillium™ / sr Integrirani neobloženi CVR/oksigenator sa  
Trillium™ Biosurface površinom / sv Integrerad obelagd CVR/oxigenator med  
Trillium™ bioyta / tr Trillium™ Biyo Yüzeyli Entegre Kaplamasız CVR/Oksijenatör

Important information prior to use: Read the instructions. Lire les instructions. Lesen Sie die Anweisungen.

Consult instructions  
for use

USA Rx Only  
For US audiences only

Store at  
controlled temperature level

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Consult instructions  
for use

USA Rx Only  
For US audiences only

Store at  
controlled temperature level

Do not use if  
package is damaged

STERILE  
Sterilized using ethylene oxide

Do not re-sterilize

Nonpyrogenic, flat path

Use

Do not reuse

Qty: 1  
Quantity

YYYY-MM-DD  
Use-by date

YYYY-MM-DD  
Date of manufacture

Barcode Area  
FOR PLACEMENT ONLY

Medtronic, Inc.  
710 Medtronic Parkway  
Minneapolis, MN 55432 USA

CE  
0123

M958674001 Rev. 1A

Affinity NT oxygenator with Balance biosurface, including  
Example of Shipper Label

## Medtronic

### Affinity NT™

REF 541B  
Catalog number

en Oxygenator with Balance™ Biosurface\* and Cardiotomy/Venous Reservoir /  
bg Оксигенатор с биоповърхност Balance™ и кардиотомия/венозен резервоар /  
cs Oxygenátor s biopovrchem Balance™ a kardiatický/venózní rezervoár /  
da Oxygenator med Balance™ biooverflade og kardiatomicreservoir/venøst reservoir /  
de Oxygenator mit Balance™ Biosurface und Kardiatomicreservoir/venöses Reservoir /  
el Οξυγονωτής με βιοεπιφάνεια Balance™ και δεξαμενή αίματος καρδιотомия/φλεβικού αίματος / es Oxigenador con superficie biocompatible Balance™ y reservorio de cardiotomia/venoso / et Balance™ biokattega oksügaator ja kardiotomia/venoosne reservuaar / fi Balance™ biopinnoitteella päällystetty hapetin ja kardiotomia-/askimoverisäiliö / fr Oxigénéateur avec biosurface Balance™ et réservoir de cardiotomie/veineux / hr Oksigenator s biopovršinom Balance™ i kardiatičkim/venskim spremnikom / hu Balance™ biofelülettel ellátott oxigénátor és kardiatiomás vagy vénás tartály / it Ossigenatore con superficie biopassiva Balance™ e cardiotomo/serbatoio venoso / lt Oksigenatorius su „Balance™“ biologine dangra ir kardiatiomijos / veninio rezervuaru / lv Oksigenators ar Balance™ bioloģisko pārklājumu un kardiatiomijas/venozais rezervuārs / mk Oksigenator co Balance™ Biosurface и кардиотомски/венски резервоар / nl Oxygenator met Balance™-bio-oppervlakte/aag en cardiotomie-/veneus reservoir / no Oksygenator med Balance™ biooverflate og kardiatomicreservoir / venøst reservoir / pl Oksygenator z biopowłoką Balance™ i zbiornik kardiatomicjny/żylny / pt-br Oxigenador com biosuperfície Balance™ e reservatório de cardiotomia/venoso / pt-pt Oxigenador com biosuperfície Balance™ e reservatório de cardiotomia/venoso / ro Oxigenator cu biosuprafață Balance™ și rezervor de cardiotomie/venos / ru Оксигенатор с биопокрытием Balance™ и кардиотомный/венозный резервуар / sk Oksigenátor s biologickou povrchovou úpravou Balance™ a kardiatickým/venóznym zásobníkom / sl Oksigenator z biološko površino Balance™ in rezervoar za kardiatiomjsko/vensko kri / sr Oksigenator sa Balance™ Biosurface površinom i kardiatiomjskim/venskim rezervoarom / sv Oxigenator med Balance™ bioyta och kardiatiomi-/venreservoir / tr Balance™ Biyo Yüzeyli Oksijenatör ve Kardiatiomi/Venöz Rezervuar / uk Оксигенатор з біопокриттям Balance™ i кардіотомний/венозний резервуар

\*Technology licensed under agreement from BioInteractions, Limited, United Kingdom.

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



Once sealed in the header bag, it will then be placed inside a vented poly bag followed by being packed into a corrugated sleeve. 4 product filled corrugated sleeves will then be packed into the shipper. 1 IFU will be included with the 4-pack configuration and placed in the shipper.

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

The combo configuration of the proposed device will be packed in a Tyvek header bag. Once sealed in the header bag will then be set in 2 molded pulp trays - one on each side of the device. The molded pulp trays will be constrained by sliding a corrugate sleeve over the product filled molded pulp trays and then secured with locking tabs. 2 product filled corrugate sleeves will then be packed into the shipper. 1 IFU will be included with the 2-pack configuration and placed in the shipper.

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

|                                                                                       |                                                                                                                                        |
|---------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
|     | Manufactured in                                                                                                                        |
|    | Serial number                                                                                                                          |
|    | Catalog number                                                                                                                         |
|    | Quantity                                                                                                                               |
|    | For US audiences only                                                                                                                  |
|    | Use By                                                                                                                                 |
|    | CE Mark. Conformité Européenne (European Conformity). This symbol means that the device fully complies with European Council Directive |
|    | Fragile                                                                                                                                |
|    | Keep dry                                                                                                                               |
|   | Upper side                                                                                                                             |
|  | Upper Limit of Temperature                                                                                                             |
|  | Transit temperature limit                                                                                                              |
|  | Attention, See instructions for Use                                                                                                    |

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

|                                                                                     |                         |
|-------------------------------------------------------------------------------------|-------------------------|
|  | Nonpyrogenic Fluid Path |
|  | Contains phthalate DEHP |

Affinity NT oxygenator with Balance biosurface, including

| Symbol                                                                                | Description                     |
|---------------------------------------------------------------------------------------|---------------------------------|
|    | Sterilized using ethylene oxide |
|    | Do not resterilize              |
|    | Lot Number                      |
|    | Transit temperature limit       |
|   | Do not reuse                    |
|  | Date of manufacture             |
|  | Manufacturer                    |
|  | Serial number                   |
|  | Quantity                        |
|  | For U.S. audiences only         |

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**3. Risk reduction “as far as possible” versus “as low as reasonably practicable”:**

- a. The team considered that the risks to be reduced as low as possible when considered in the context of the highly skilled clinical perfusionists using the device during cardiopulmonary bypass procedures.

**4. Discretion as to whether a risk-benefit analysis needs to take place:**

- a. The team assessed the overall risks associated with the device and determined that the benefits of using the device outweigh the risks as documented in both the DFMECA and the Risk Management Report.

**5. Discretion as to the risk control options/measures:**

- a. Assessment and application of the risk mitigations was completed as outlined in Section 6.2 of EN ISO 14971.

**6. Deviation as to the first risk control option:**

- a. The risk control options were applied as outlined in Section 6.2 of EN ISO 14971 with the emphasis placed on designing out risks.

**7. Information of the users influencing the residual risk:**

- a. The Instructions for Use (IFU) were not used as the sole risk reduction measure during the risk analysis process.
- b. The IFU is used only to inform the user of the residual risks related to the device.

Figure 4 represents the Risk Evaluation/Risk Classification Matrix applied for all Structural Heart devices. A zone of risk is assigned based on the point where the occurrences intersect on the matrix.

Risks are evaluated in the following manner, depending on the zone:

|        |                                                                                                                                                                                                                                                         |
|--------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Zone 3 | Reduce risks within this region to a Zone 2 or Zone 1 risk with Risk Control measures. Is a risk reduction is not feasible, document that the risk has been reduced to as low as possible and evaluate whether the benefit outweighs the residual risk. |
|--------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

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
Reduce all Zone 2 risks with Risk Control measures. If the Zone 2 risk directly results from the design and/or manufacturing failures, pursue additional controls or design improvements to reduce the risk. If none of these approaches are feasible, document that the risk has been reduced as low as possible and evaluate whether the benefit outweighs the residual risk.

Within this region the level of risk may be considered acceptable once the risk has been reduced as low as possible. At a minimum, document existing risk controls.

**RISK EVALUATION / RISK CLASSIFICATION MATRIX**


**Figure 4: Risk Evaluation Matrix**

The device performance risks (functional hazards) and material-related risks (for example, biocompatibility) that were assessed are independent of the intended use by clinician (the functional risks of the SSA and NSB versions are the same).

#### **Affinity NT oxygenator with Balance biosurface**

Each device was assessed in accordance with EN ISO 14971:2009. In addition, a gap analysis was conducted for compliance to EN ISO 14971:2012 and because the technical contents of the standard are unchanged from EN ISO 14972:2009 compliance to the updated standard is considered to be maintained. The seven content

**RISK EVALUATION / RISK CLASSIFICATION MATRIX**  
**OCCURRENCE RANK**

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devia within 14971 Annex ZA (informative) outlining the relationship between EN ISO 14971:2012 and the requirements of the EU Directive 93/42/EEC on medical devices are addressed as follows:

1. Treatment of negligible risks:

a. Review of the DFMECA background section of each device documented that the team chose to mitigate all risks.

b. Review of the overall residual risks demonstrated a positive risk benefit profile when considering the devices are used in clinical practice by highly skilled clinical perfusionists who are thoroughly trained in cardiopulmonary bypass procedures.

2. Discretionary power of manufacturers as to the acceptability of risks:

a. The post market surveillance information is reviewed and updated on an ongoing basis including quarterly reviews to ensure the risks are known and properly documented.

The post market surveillance report for oxygenator systems (10007385DOC) documents the ongoing acceptability of risk.

b. The team considered the acceptability of the risk to be equivalent to the devices currently on the market.

c. All risks have been taken into account when considering the overall risk benefit of the oxygenators and oxygenator systems.

3. Risk reduction "as far as possible" versus "as low as reasonably practicable":

a. The team considered that the risks to be reduced as low as possible when considered in the context of the highly skilled clinical perfusionists using the device during cardiopulmonary bypass procedures. Economic considerations were not part of the risk control measures.

4. Discretion as to whether a risk-benefit analysis needs to take place:

a. The team assessed the overall risks associated with the device and determined that the benefits of using the device outweigh the risks as documented in both the DFMECA and the Risk Management Report.

5. Discretion as to the risk control options/measures:

a. Assessment and application of the risk mitigations was completed as outlined in Section 6.2 of EN ISO 14971.

b. Emphasis was placed on designing out the risks and effectiveness of these design controls was reviewed and assessed throughout all phases of the validation process. The verification and validation activities are summarized in Section 3 of this technical document that cover pre-clinical testing.

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8 The conclusion of this qualification report is that optimized Trillium coating

process for manufacturing Affinity NT Hollow Fiber Oxygenators is considered qualified based on successfully meeting all test requirement.

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1 Table I below is a summary of the qualification test results for gas transfer characteristics.

| TABLE I<br>Summary of Gas<br>Transfer<br>Characteristics<br>(Q <sub>blood</sub> = 7.0 LPM) |                   |                                      |                                       |                              |
|--------------------------------------------------------------------------------------------|-------------------|--------------------------------------|---------------------------------------|------------------------------|
| Sample I.D.                                                                                | # Units<br>Tested | O <sub>2</sub> Transfer<br>(ml/min.) | CO <sub>2</sub> Transfer<br>(ml/min.) | Pressure Drop<br>(mmHg)      |
| Standard Trillium                                                                          | 22                | Mean = 387<br>StDev = 8.0            | Mean = 305.5<br>StDev = 12.4          | Mean = 103.8<br>StDev = 6.9  |
| HH (1/5 conc.)                                                                             | 22                | Mean = 398.7<br>StDev = 8.7          | Mean = 316.9<br>StDev = 15.7          | Mean = 105.8<br>StDev = 10.8 |
| LL (1/7 conc.)                                                                             | 22                | Mean = 399.1<br>StDev = 8.9          | Mean = 322.1<br>StDev = 8.7           | Mean = 105.6<br>StDev = 7.1  |
| product specification                                                                      |                   | ≥ 374 ml/min.                        | ≥ 288 ml/min.                         | < 124 mmHg                   |

2 Medstat 4.3- Oneway-ANOVA statistical model was used to test the null hypothesis that the test (1/5 and 1/7 cone.) means and standard deviations vs. standard Trillium control means and standard deviations are equivalent at the 95% confidence interval. The model also provides for an alternative hypothesis that at least two of the means are significantly different compared to the other populations.

2.1 The test results of the optimized Trillium coated units (1/5 and 1/7 cone.) show a positive significant difference for oxygen transfer and carbon dioxide transfer compared to standard Trillium coated control units. This positive difference is evident in the higher mean gas transfer levels for the optimized Trillium coated units which translates into overall better gas transfer performance and therefore passes the defined acceptance criteria of equivalency (or better) to control units.

2.2 No statistical differences were detected between optimized Trillium coated test and control units for pressure drop across the Oxygenator.

3 Table II below is a summary of the Blood Trauma test results reported as a percent of the initial corrected value.

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3.1 Comparing the mean percents of initial value for Blood Trauma of the optimized Trillium coated units (1/5 & 1/7 concentration) vs. standard Trillium coated units shows functionally equivalent test results were obtained throughout the six (6) hour test duration.

4 Table III below is a summary of the qualification test results for Plasma Breakthrough.

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| TABLE III<br>Summary of Plasma Breakthrough |         |                                               |               |
|---------------------------------------------|---------|-----------------------------------------------|---------------|
| Sample I.D.                                 | # Units | Breakthrough Time (min.)                      | Accept Result |
| Uncoated Controls                           | 7       | Mean = 139.1                                  | Pass          |
| Standard Trillium                           | 7       | No Breakthrough                               | Pass          |
| HHH (1/5 conc.)                             | 8       | (n=6); No Breakthrough<br>(n=2); Mean = 337.5 | Pass          |
| LLL (1/7 conc.)                             | 6       | (n=4); No Breakthrough<br>(n=2); Mean = 345   | Pass          |
| product specification                       |         | ≥ 30min.                                      |               |

4.1 The Plasma breakthrough test results show comparable breakthrough levels were obtained for the HHH (1/5) and LLL 1/7) test units to standard Trillium test units. All seven (7) of the standard Trillium test units exhibited no breakthrough compared to six (6) of the eight (8) units for the HHH and four (4) of the six (6) units exhibited no breakthrough. The remaining two (2) units each of HHH and LLL breakthrough times were 337.5 and 345 minutes respectively.

5 Table IV below is a summary of the qualification test results for Leaching.

Gel Permeation Chromatograms for the individual test units indicate that no leachables are coming out of the materials from the HHH (1/5) or LLL (1/7) coated units as evidence from the GPC results of no peaks detected in the range of 20 — 45 minutes.

| TABLE IV<br>Summary of Leach Results |         |                   |             |
|--------------------------------------|---------|-------------------|-------------|
| Sample I.D.                          | # Units | GPC Results       | Test Result |
| Uncoated Control                     | 5       | No Peaks Detected | PASS        |
| HHH (1/5 conc.)                      | 5       | No Peaks Detected | PASS        |
| LLL (1/7 conc.)                      | 5       | No Peaks Detected | PASS        |

6 Table V below is a summary of the qualification test results for Biocompatibility (includes cytotoxicity, hemolysis, pyrogen and bioburden testing). The test results indicate no biocompatibility issues are evident with either of the HHH (1/5 conc.) or LLL (1/7 conc.) coated units.

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| TABLE V                             |         |              |               |                  |                |
|-------------------------------------|---------|--------------|---------------|------------------|----------------|
| Summary of Biocompatibility Testing |         |              |               |                  |                |
| Sample I.D.                         | # Units | Cytotoxicity | Hemolysis     | Bioburden        | Pyrogen        |
| HHH (1/5 conc.)                     | 5       | Grade 0      | Non-hemolytic | < 18             | < 1.4          |
| Test Result                         |         | PASS         | PASS          | PASS             | PASS           |
| LLL (1/7 Conc.)                     | 5       | Grade 0      | Non-hemolytic | < 18             | < 1.4          |
| Test Result                         |         | PASS         | PASS          | PASS             | PASS           |
| product specification               |         | ≤ Grade 2    | Non-hemolytic | < 400 CFU/device | < 20 EU/device |

## CONCLUSION

The conclusion of this qualification report is that optimized Trillium coating process for manufacturing Affinity NT Hollow Fiber Oxygenators is considered qualified based on successfully meeting all test requirement.

Report for the Validation of Optimized Trillium Coating Process of Affinity NT Hollow Fiber Oxygenator

## SUMMARY OF RESULTS

1 The optimized validation testing was acceptable for all tests performed. The optimized coating process is considered qualified for use in production.

2 The Trillium coated Affinity NT Hollow Fiber Oxygenators produced with the optimized process, (model 511T) were produced in accordance with test protocol. The validation was done at two concentrations of polymer concentrations, 1/4 the current concentration for non-aged testing and 1/7 as a worse case for shelf life testing.

3 The Optimized Trillium coated Affinity NT hollow fiber oxygenators were coated per manufacturing procedures 31410 (1/4 polymer cone.) and 31383 (1/7 polymer cone.) respectively. The standard Trillium coated Oxygenators were coated per manufacturing procedure 31114; Affinity HFO, Trillium.

4 The validation test results demonstrated that the Trillium coated oxygenators coated with the optimized process showed consistency for gas transfer characteristics meeting the criteria in product specification 151034-01 and the defined acceptance criteria of equivalency (or better) to the standard Trillium control group. Gas transfer test results show increased gas transfer that is similar to the increase that has been seen in the base fiber used in the time frame of the control units to the fiber used for the test units. Gas transfer characteristics include oxygen transfer, carbon dioxide transfer and blood pressure drop across the coated oxygenator. Refer to Table I - Summary of Gas Transfer Characteristics.

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5 The validation test results demonstrated that the Trillium coated Affinity oxygenators coated with the optimized process showed functional consistency between the lots and yielded acceptable product that shows functionally equivalent blood trauma results to the standard Trillium coated and uncoated control units. Refer to Tables II - IV - Summary of Blood Trauma.

6 The validation test results demonstrated that the Trillium coated Affinity oxygenators coated with the optimized process showed consistency between the lots and yielded acceptable product for plasma breakthrough meeting the criteria in product specification 151034-01 of > 30 min. time to breakthrough fiber. Refer to Table V - Summary of Plasma Breakthrough.

7 The validation testing of the oxygenators coated with the optimized process resulted in no leaching in either the aged units (at 1/7 concentration) or non-aged (at concentration) evident that is caused by the presence of the coating materials. Refer to Table VI and VIII - Summary of Leaching Data.

8 The validation test of the dye coverage of the oxygenators coated with the optimized process showed complete coverage for all aged and non-aged units.

See results in Table VII and IX - Summary of dye coverage results.

9 The validation test results demonstrated that the Trillium oxygenators coated with the optimized process and accelerated aged for 18 months yielded acceptable product for biocompatibility characteristics per individual Geneva Lab test results and therefore meeting the defined acceptance of the test plan. Biocompatibility characteristics include cytotoxicity and hemolysis.

#### QUALIFICATION/VALIDATION RESULTS

1 Table I below is a summary of the validation test results for gas transfer characteristics.

| TABLE I<br>Summary of Gas<br>Transfer<br>Characteristics<br>(Q <sub>blood</sub> = 7.0 LPM) |                   |                                      |                                       |                         |
|--------------------------------------------------------------------------------------------|-------------------|--------------------------------------|---------------------------------------|-------------------------|
| Sample I.D.                                                                                | # Units<br>Tested | O <sub>2</sub> Transfer<br>(ml/min.) | CO <sub>2</sub> Transfer<br>(ml/min.) | Pressure Drop<br>(mmHg) |
| Standard Trillium                                                                          | 22                | 378.52<br>9.9                        | 296.18<br>14.1                        | 102.32<br>5.8           |
| Batch 1                                                                                    | 22                | 395.5<br>6.6                         | 310.32<br>11.4                        | 108.73<br>6.0           |
| Batch 2                                                                                    | 22                | 392.22<br>7.0                        | 304.36<br>10.9                        | 113.23<br>5.6           |
| Batch 3                                                                                    | 22                | 390.08<br>7.4                        | 305.89<br>12.8                        | 110.32<br>7.7           |
| Product specification                                                                      |                   | ≥ 374 ml/min.                        | ≥ 288 ml/min.                         | < 124 mmHg              |

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2 All units met the specification condition for oxygen transfer, carbon dioxide transfer and pressure drop.

3 Statistical analysis (from Medstat 4.3) comparing the three validation lots against the control group. The analysis shows a higher oxygen transfer, carbon dioxide transfer and pressure drop between the test and control groups and a difference in oxygen transfer between validation runs 1 and 3 at the 95% confidence interval.

3.1 The difference between the test and control groups gas transfer results is positive for oxygen transfer and carbon dioxide transfer and negative for pressure drop. The control group oxygenators were taken from finished goods (built in January and February). The base performance of the fiber has increased from this time to the current production (as well as increasing the pressure drop). The increase in performance from the control units (older fiber) to the test units (newer fiber) can be accredited to the changes in the base fiber performance. The important fact is that all units passed all specification criteria.

4 Table II below is a summary of the Blood Trauma test results reported as a percent of the initial corrected value.

| Table II: Percent White Blood Cell reduction |                 |               |               |               |                    |               |               |               |
|----------------------------------------------|-----------------|---------------|---------------|---------------|--------------------|---------------|---------------|---------------|
| Time (hr)                                    | Average         |               |               |               | Standard deviation |               |               |               |
|                                              | Tri1<br>control | Val.<br>Run 1 | Val.<br>Run 2 | Val.<br>Run 3 | Tri1<br>control    | Val.<br>Run 1 | Val.<br>Run 2 | Val.<br>Run 3 |
| 0,5                                          | 79.83           | 74.67         | 73.33         | 74.67         | 11.16              | 9.07          | 13.26         | 11.83         |
| 1                                            | 75.17           | 72            | 70.67         | 71.83         | 14.13              | 8.367         | 10.42         | 8.796         |
| 2                                            | 66.5            | 66            | 65.33         | 66.67         | 7.868              | 5.441         | 9.75          | 7.528         |
| 3                                            | 64.17           | 58            | 57            | 61.33         | 13.29              | 9.757         | 10.06         | 8.847         |
| 4                                            | 61.67           | 57            | 59.67         | 59            | 8.959              | 5.215         | 9.873         | 8.39          |
| 5                                            | 62.83           | 59.33         | 64.33         | 56.67         | 12.58              | 5.391         | 13.71         | 4.082         |
| 6                                            | 59.5            | 57.33         | 60.33         | 60.17         | 9.247              | 7.062         | 6.89          | 8.208         |

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4.1 Comparing the mean change for the validation test units vs. standard Trillium coated units shows functionally equivalent test results.

5 Table V below is a summary of the validation test results for Plasma Breakthrough.

| TABLE V<br>Summary of Plasma Breakthrough |         |                          |               |
|-------------------------------------------|---------|--------------------------|---------------|
| Sample I.D.                               | # Units | Breakthrough Time (min.) | Accept Result |
| Uncoated Controls                         | 7       | (n=5); Mean = 222.2      | Pass          |
| Standard Trillium                         | 7       | (n=1); No Breakthrough   | Pass          |
| Batch 1                                   | 7       | (n=7); No Breakthrough   | Pass          |
| Batch 2                                   | 7       | (n=7); No Breakthrough   | Pass          |
| Batch 3                                   | 7       | (n=1); Mean = 270        | Pass          |
|                                           |         | (n=6); No Breakthrough   | Pass          |
| Product specification                     |         | ≥ 30min.                 |               |

5.1 All units (test, Trillium control and uncoated control) had significantly longer plasma breakthrough times that the specification time of 30 minutes. A qualitative inspection shows that the uncoated control units had a shorter plasma breakthrough time (based on 6 of the 7 units with breakthrough during the test time versus only 1 of the 21 validation units).

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6 Table VI below is a summary of the validation test results for leaching test. Gel Permeation Chromatograms for the individual test units indicate that no leachables are coming out of the materials from the coated units as evidence from the GPC results of no peaks detected in the range of 20 - 45 minutes.

| Table VI: Summary of Leach Results (1/4 concentration from TP-0077-003) |            |                   |             |
|-------------------------------------------------------------------------|------------|-------------------|-------------|
| Test group                                                              | # of units | GPC Results       | Test Result |
| Uncoated control                                                        | 5          | No peaks detected | Pass        |
| Validation run 1                                                        | 5          | No peaks detected | Pass        |
| Validation run 2                                                        | 5          | No peaks detected | Pass        |
| Validation run 3                                                        | 5          | No peaks detected | Pass        |

7. Table VII below is a summary of the validation test results for the dye coverage test for non-aged units at 1/4 current concentration. Inspection of the units after dye test shows uniform complete coverage on all test units.

| Table VII: Summary of Dye test Results (1/4 concentration from TP-0077-003) |            |                      |             |
|-----------------------------------------------------------------------------|------------|----------------------|-------------|
| Test group                                                                  | # of units | Dye coverage results | Test Result |
| Uncoated control                                                            | 5          | Uniform coverage     | Pass        |
| Validation run 1                                                            | 5          | Uniform coverage     | Pass        |
| Validation run 2                                                            | 5          | Uniform coverage     | Pass        |
| Validation run 3                                                            | 5          | Uniform coverage     | Pass        |

8 Table VIII below is a summary of the validation test results for aged unit leaching test at 1/7 current concentration. Gel Permeation Chromatograms for the individual test units are provided in Attachment F and indicate that no leachables are coming out of the materials from the coated units as evidence from the GPC results of no peaks detected in the range of 20 - 45 minutes.

| Table VIII: Summary of Leach Results (aged units at 1/7 concentration from TP-0043-124) |            |                   |             |
|-----------------------------------------------------------------------------------------|------------|-------------------|-------------|
| Test group                                                                              | # of units | GPC Results       | Test Result |
| Validation run 1                                                                        | 4          | No peaks detected | Pass        |
| Validation run 2                                                                        | 7          | No peaks detected | Pass        |
| Validation run 3                                                                        | 4          | No peaks detected | Pass        |

9 Table IX below is a summary of the validation test results for the dye coverage

test for aged units at 1/7 current concentration. Inspection of the units after dye test shows uniform complete coverage on all test units.

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| Table IX: Summary of Dye Coverage Results (aged units at 1/7 concentration from TP-0043-124) |            |                      |             |
|----------------------------------------------------------------------------------------------|------------|----------------------|-------------|
| Test group                                                                                   | # of units | Dye coverage Results | Test Result |
| Validation run 1                                                                             | 4          | Uniform coverage     | Pass        |
| Validation run 2                                                                             | 7          | Uniform coverage     | Pass        |
| Validation run 3                                                                             | 4          | Uniform coverage     | Pass        |

10 Table X below is a summary of the validation test results for aged unit Biocompatibility (includes cytotoxicity and hemolysis testing). The test results indicate no biocompatibility issues are evident with the validation coated units.

11 The sample size of the hemolysis testing was supposed to be 5 units per lot. Because of integrity failures that are unrelated to the Trillium process, the number of samples available for the test from validation run 1 was only 1 unit, and for validation run 2 was only 4 units. The total sample size was obtained using additional units from validation run 3. This substitution is acceptable since the definition of a validation run is 60 consecutive coated units. The units are actually produced as a 10 unit run, so a validation run is 6 consecutive runs. The validation build is 18 coating runs with the first 6 called validation run 1 and so on.

| Table X: Summary of Biocompatibility Testing (aged units from TP-0043-124) |              |              |            |              |               |            |
|----------------------------------------------------------------------------|--------------|--------------|------------|--------------|---------------|------------|
| Run number                                                                 | No. of units | Cytotoxicity |            | Hemolysis    |               |            |
|                                                                            |              | Test result  | Acceptance | No. of units | Test result   | Acceptance |
| Validation run 1                                                           | 5            | Grade 0      | Pass       | 1            | Non-hemolytic | Pass       |
| Validation run 2                                                           | 5            | Grade 0      | Pass       | 4            | Non-hemolytic | Pass       |
| Validation run 3                                                           | 5            | Grade 0      | Pass       | 10           | Non-hemolytic | Pass       |

## CONCLUSION

The conclusion of this validation report is that optimized Trillium coating process for manufacturing Affinity NT Hollow Fiber Oxygenators is considered validated based on successfully meeting all test requirements as detailed in test protocol.

Validation of Anesthesia Use (Isoflourane, Sevoflourane) with the Affinity HFO, Report

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#### SUMMARY OF RESULTS

This report contains the test results of the Trillium coated Affinity HFO (511T).

1 The results of the isoflurane testing passed the requirements.

2 The results of the sevoflurane testing passed the requirements.

#### REQUIREMENTS AND ACCEPTANCE CRITERIA

1 Oxygen Transfer Rate

The oxygen transfer rate shall be a mean of  $> 374$  ml  $O_2$ /min. at  $7 \pm .3$  l/min blood flow, and a gas flow rate of  $7 \pm .5$  l/min., under venous inlet conditions (Table 1), Test Method #31429

2.2 Carbon Dioxide Transfer Rate

The carbon dioxide transfer rate shall be a mean of  $> 288$  ml  $CO_2$ /min at  $7 \pm 0.3$  L/min

blood flow, and a gas flow rate of  $7 \pm .5$  L/min under venous inlet conditions (Table 1), Test Method #31429

Table 1: Venous Inlet Conditions

|                          |                    |
|--------------------------|--------------------|
| Blood Temperature        | $37 \pm 2^\circ C$ |
| $O_2$ Saturation         | $65 \pm 2\%$       |
| $pCO_2$                  | $45 \pm 2$ mmHg    |
| Hemoglobin Concentration | $12 \pm .3$ gm/dl  |
| Base Excess              | $0 \pm .5$ meq/l   |

3 Mean Blood Pressure Drop

The mean blood pressure drop across the Oxygenator shall not exceed 124 mmHg at a blood flow rate of  $7.0 \pm .3$  L/min, Test Method #31429

4 Mean Gas Pressure Drop

The mean gas pressure drop across the Oxygenator shall not exceed 15 mmHg at a ventilating gas flow rate of  $7.0 \pm .5$  L/min, Test Method #31429

5 Pressure Integrity

The Oxygenator assembly shall not produce a leak when subjected to mechanical integrity testing in accordance with Test Method #31220

6 Plasma Breakthrough

The relative plasma breakthrough time shall be 30 minutes minimum, Test Method #31432

Note: All tolerance interval calculations will be at a 95% confidence level with a 90% reliability level

#### VALIDATION RESULTS

The following tables summarize the test result data found in Attachments 1-3\*. This data was then analyzed using Medstat v5.2 to determine the descriptive statistics and the 95% / 90% confidence intervals where applicable\*\*.

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\*Note: The gas transfer rate and the pressure drop tests (sections 6.1-6.4) contain data from "time dependant" (sample# 1-3) and "sweeps" testing (sample# 4-6).

\*\*Note: The "Oxygen Transfer Rate" and "Carbon Dioxide Transfer Rate" performance acceptance criteria are based on the means being above a minimum value; therefore, the use of confidence intervals is not applicable.

#### 1 Oxygen Transfer Rate

Table 2: Oxygen Transfer Rate Results (ml O<sub>2</sub>/min)

| Sample # | Isoflurane |        | Sevoflurane |        |
|----------|------------|--------|-------------|--------|
|          | 0 hrs      | 6 hrs  | 0 hrs       | 6 hrs  |
| Control  | 396.95     | 375.19 | 411.87      | 390.27 |
| 1        | 385.81     | 369.04 | 396.98      | 397.64 |
| 2        | 389.74     | 372.83 | 407.27      | 401.45 |
| 3        | 384.90     | 335.83 | 402.02      | 388.47 |
| 4        | 364.44     | -      | 369.87      | -      |
| 5        | 381.55     | -      | 358.11      | -      |
| 6        | 375.20     | -      | 366.07      | -      |
| Average  | 380.27     | 359.23 | 383.39      | 395.85 |

#### 2 Carbon Dioxide Transfer Rate

Table 3: Carbon Dioxide Transfer Rate Results (ml CO<sub>2</sub>/min)

| Sample # | Isoflurane |        | Sevoflurane |        |
|----------|------------|--------|-------------|--------|
|          | 0 hrs      | 6 hrs  | 0 hrs       | 6 hrs  |
| Control  | 298.43     | 266.84 | 292.71      | 345.83 |
| 1        | 310.24     | 272.42 | 295.77      | 292.35 |
| 2        | 301.91     | 255.36 | 311.16      | 290.98 |
| 3        | 303.98     | 226.70 | 302.77      | 286.20 |
| 4        | 296.02     | -      | 298.93      | -      |
| 5        | 292.02     | -      | 281.67      | -      |
| 6        | 306.85     | -      | 283.89      | -      |
| Average  | 301.84     | 251.49 | 295.70      | 289.84 |

#### 3 Mean Blood Pressure Drop

Table 4: Mean Blood Pressure Drop Results (mmHg)

| Sample # | Isoflurane    |           | Sevoflurane   |           |
|----------|---------------|-----------|---------------|-----------|
|          | Pressure Drop | Test Date | Pressure Drop | Test Date |
| 1        | 104           | 2-28-08   | 81            | 2-26-08   |
| 2        | 100           |           | 95            |           |
| 3        | 106           |           | 105           |           |
| 4        | 117           | 3-5-08    | 116           | 3-5-08    |
| 5        | 117           |           | 117           |           |
| 6        | 120           |           | 115           |           |

The "mean blood pressure drop" tests were completed over several days and blood varies from day-to-day as evidenced by table 4.

The primary variation for the blood pressure drop is the blood, not the presence or

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absent anesthesia. Past quarterly testing of Trillium coated product (Q4CY07, 11-13-07) shows that the blood pressure drop (96 to 104 mmHg) is approximately the same as the test data obtained for anesthesia testing.

Due to the fact the primary source of blood pressure drop comes from the blood, the "mean blood pressure drop" acceptance criteria is not applicable for anesthesia evaluation.

Since blood varies from day-to-day, the sample size is reduced to down to three (n=3), which is too small for effective statistical analyses; therefore, confidence intervals will not be calculated. Note: all samples did meet the blood pressure drop specification of  $\leq 124$  mmHg.

#### 4 Mean Gas Pressure Drop

Table 5: Mean Gas Pressure Drop Results (mmHg)

| Sample # | Isoflurane | Sevoflurane |
|----------|------------|-------------|
| 1        | 7.7        | 8.1         |
| 2        | 8.0        | 7.8         |
| 3        | 7.4        | 7.5         |
| 4        | 7.9        | 8.0         |
| 5        | 8.3        | 7.7         |
| 6        | 8.2        | 7.7         |

Table 6: Mean Gas Pressure Drop Descriptive Statistics

| Descriptive Statistics | Isoflurane | Sevoflurane |
|------------------------|------------|-------------|
| N                      | 6          | 6           |
| Mean                   | 7.91667    | 7.80000     |
| Std Dev                | 0.33116    | 0.21909     |
| Variance               | 0.10967    | 0.04800     |
| Range                  | 0.90000    | 0.60000     |
| Skewness               | -0.56539   | 0.17116     |
| Kurtosis               | -0.33952   | -0.78125    |
| Minimum                | 7.40000    | 7.50000     |
| 25% Point              | 7.70000    | 7.70000     |
| Median                 | 7.95000    | 7.75000     |
| 75% Point              | 8.20000    | 8.00000     |
| Maximum                | 8.30000    | 8.10000     |
| Test for Normality     |            |             |
| S-W                    | 0.96611    | 0.95517     |
| P-Value                | 0.86534    | 0.79538     |

Do not reject normality assumption at 0.05 level significance  
Warning: Sample size too small to accept normality.

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Table 7: Mean Gas Pressure Drop Tolerance Interval Calculations

|                                                          |            |             |
|----------------------------------------------------------|------------|-------------|
| Tolerance Interval<br>Normal - Calculate Interval Limits |            |             |
| Inputs                                                   | Isoflurane | Sevoflurane |
| Confidence (%)                                           | 95         | 95          |
| % of Population Included in Interval                     | 90         | 90          |
| Sample Size                                              | 6          | 6           |
| Mean                                                     | 7.91667    | 7.80000     |
| Standard Deviation                                       | 0.33116    | 0.21909     |
| Sided                                                    | Upper      | Upper       |
| Output                                                   |            |             |
| Upper Tolerance Limit                                    | 8.91222    | 8.45864     |

5 Pressure Integrity

Table 8: Pressure Integrity Results

| Isoflurane    | Sevoflurane   |
|---------------|---------------|
| 29 of 29 PASS | 29 of 29 Pass |

6 Plasma Breakthrough

Table 9: Plasma Breakthrough Results

| Sample # | Isoflurane | Sevoflurane |
|----------|------------|-------------|
| 1        | n/a*       | n/a*        |
| 2        | n/a*       | n/a*        |
| 3        | n/a*       | n/a*        |
| 4        | n/a*       | n/a*        |
| 5        | n/a*       | n/a*        |

\*Note: No breakthrough after 6hrs.

Since there were no plasma breakthrough failures, descriptive statistics and tolerance interval calculations can not be determined.

## CONCLUSION

The following table summarizes the test results as it relates to the performance specification:

Table 10: Results Summary

| Test Description             | Requirement                            | Isoflurane                               | Sevoflurane                              |
|------------------------------|----------------------------------------|------------------------------------------|------------------------------------------|
| Oxygen Transfer Rate         | Ave $\geq 374$ ml O <sub>2</sub> /min  | PASS<br>Ave = 380ml O <sub>2</sub> /min  | PASS<br>Ave = 393ml O <sub>2</sub> /min  |
| Carbon Dioxide Transfer Rate | Ave $\geq 288$ ml CO <sub>2</sub> /min | PASS<br>Ave = 302ml CO <sub>2</sub> /min | PASS<br>Ave = 296ml CO <sub>2</sub> /min |
| Mean Blood Pressure Drop     | 95%/90% $\leq 124$ mmHg                | Not Performed                            | Not Performed                            |
| Mean Gas Pressure Drop       | 95%/90% $\leq 15$ mmHg                 | PASS<br>8.9 mmHg                         | PASS<br>8.5 mmHg                         |
| Pressure Integrity           | 29 of 29 PASS                          | PASS<br>29 of 29                         | PASS<br>29 of 29                         |
| Plasma Breakthrough          | 95%/90% $\geq 30$ min                  | PASS<br>$>360$ min**                     | PASS<br>$>360$ min**                     |

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\*Note: The "mean blood pressure drop" is not an applicable test for anesthesia testing, as previously discussed.

\*\*Note: No breakthrough after 6hrs.

The Trillium coated Affinity HFO (511T) samples were able to meet the performance specifications for gas transfer, pressure integrity, and plasma break through when used with Isoflurane and Sevoflurane. The time dependant gas transfer rates for the isoflurane and sevoflurane samples were about the same as the control samples; the presence of anesthesia during testing does not appear to degrade the performance beyond the normal wear of the device.

This study shows that the Trillium Coated Affinity NT HFO can be indicated for use with Isoflurane and Sevoflurane.

Test Study of Iso/Sevoflurane Delivery and Safety through Affinity Trillium Hollow Fiber Oxygenator

#### PURPOSE/OBJECTIVES

This test study summarizes test plans and test reports to assess the effectiveness and safety of Trillium-coated Affinity Oxygenator as a delivery device for anesthetic gas.

#### CONCLUSION

#### Conclusion of Evaluation of Isoflurane and Sevoflurane in Trillium Affinity Oxygenators

Isoflurane and sevoflurane delivered via Trillium-coated Affinity oxygenator can be detected in the blood at similar levels as when delivered by inhalation. Delivery via oxygenator results in adequate depth of anesthesia as assessed by a noxious stimulus test.

#### Conclusion of Toxicological Evaluation of Isoflurane and Sevoflurane in Trillium Affinity Oxygenators (In vitro/Non-blood Environment) Biological Evaluation

There is no evidence to indicate that the use of these two anesthetic gases in the Trillium Affinity Oxygenator at their intended concentrations would result in the formation of a new toxic chemical. In addition, negative cytotoxicity testing results, i.e. non-cytotoxic; demonstrate that no toxic species were produced when these two anesthetic gases were used. Based on these findings, the use of these anesthetic gases in this oxygenator would not present a health risk to a patient from the formation of a toxic chemical species.

#### Biocompatibility

| Description   | Standard     | Result |
|---------------|--------------|--------|
| Cytotoxicity  | ISO 10993-5  | PASS   |
| Sensitization | ISO 10993-10 | PASS   |

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|                                      |              |      |
|--------------------------------------|--------------|------|
| Irritation / Intracutaneous Toxicity | ISO 10993-10 | PASS |
| Systemic toxicity                    | ISO 10993-11 | PASS |
| Hemocompatibility                    | ISO 10993-4  | PASS |

#### Affinity NT oxygenator with Balance biosurface

##### Performance Testing - Bench

##### Background

Medtronic develops products according to a Design Control procedure. The procedure defines how a product design is controlled through the development cycle and ensures compliance with ISO/QSR requirements.

##### Sample Conditioning

Testing was completed on the following devices to show that they maintained performance over their designated shelf life:

- Unaged (t = 0)
- Accelerated aged (t = 2 YR AA)

Along with the test samples, predicate control devices (wherever required) were also subject to an identical set of conditioning steps. By subjecting the predicate control devices to the same conditioning as the test devices, the design verification was able to test for substantially equivalent performance of the test device to the predicate when subjected to the same conditions.

See Table 18-2 for a detailed list of conditioning:

Table 18-2: Sample Conditioning

| t = 0 (Unaged)                                       | t = 2 YR AA                                          |
|------------------------------------------------------|------------------------------------------------------|
| Build                                                | Build                                                |
| EtO Sterilization (2x Balance)                       | EtO Sterilization (2x Balance)                       |
| Environmental Condition (Temp. Cycle)                | Environmental Condition (Temp. Cycle)                |
| ASTM D4169 Simulated Distribution Testing            | ASTM D4169 Simulated Distribution Testing            |
| Isoflurane and Sevoflurane Conditioning <sup>1</sup> | Two Year Simulated Aging (65 days @ 55°C)            |
| Test                                                 | Isoflurane and Sevoflurane Conditioning <sup>1</sup> |
|                                                      | Test                                                 |

<sup>1</sup> Isoflurane and Sevoflurane conditioning was only done for select tests.

#### EtO Sterilization

All DV samples were 2x ethylene oxide (EtO) sterilized per Medplast cycle 604 Rev 3. This cycle takes "worst case" sterilization parameters into consideration and thus incorporates worstcase tolerances for temperature, pressure, rate of pressure change, humidity, EO concentration, and exposure across all applicable production sterilization cycles.

#### Environmental Conditioning

- Two Years Simulated Aging

Applicable samples were accelerated aged for a simulated total of two years (730 days); q = 2.0, an accelerated aging temperature of

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$55^{\circ}\text{C}$ , a product storage temperature of  $20^{\circ}\text{C}$  were used to calculate an accelerated exposure time of 65 days. This applies to the two-year accelerated aged test population.

$$T = t \cdot q^{\left(\frac{T_{AA} - T_{RT}}{10}\right)}$$

Where:

$T = 730$  days = Required shelf life

$t$  = Accelerated exposure time

$q = 2$  = Acceleration factor (typically 2.04)

$T_{AA} = 55^{\circ}\text{C}$  = Accelerated aging temperature in  $^{\circ}\text{C}$

$T_{RT} = 20^{\circ}\text{C}$  = Product storage temperature in  $^{\circ}\text{C}$

- Temperature Cycling

Samples also underwent temperature conditioning to expose finished product to simulated extreme but realistic shipping/storage temperatures.

These conditioning steps are outlined in Table 18-3.

Table 18-3: Temperature Cycling Parameters

| Sequence | Condition | Temperature Set Point       | % Relative Humidity | Exposure Time (Allowed Variability)                              |
|----------|-----------|-----------------------------|---------------------|------------------------------------------------------------------|
| 0        | Ambient   | $23 \pm 3^{\circ}\text{C}$  | $50 \pm 10\%$       | $\geq 24$ hours                                                  |
| 1        | Freeze    | $-18 \pm 2^{\circ}\text{C}$ | Uncontrolled        | 24 hours, ( $\pm 2$ hours)                                       |
| 2        | Tropical  | $40 \pm 2^{\circ}\text{C}$  | $90 \pm 5\%$        | 24 hours, ( $\pm 2$ hours)                                       |
| 3        | Desert    | $55 \pm 2^{\circ}\text{C}$  | $25 \pm 5\%$        | 24 hours, ( $\pm 2$ hours)                                       |
| 4        | Ambient   | $23 \pm 3^{\circ}\text{C}$  | $50 \pm 10\%$       | NA, sequence can end once final condition is achieved. No limit. |

Ramp rates between exposures shall be  $5^{\circ}\text{C} (\pm 1^{\circ}\text{C})$  per hour.

Distribution Conditioning

Applicable samples were shake-and-vibration conditioned per ASTM D4169-16: DC13, Assurance Level 1. This cycle is shown in Table 18-4 below with the following exclusions.

- The packaging and the product are breathable; thus, low pressure testing Schedule I was omitted
- The shipper box is not under 275 lbf/in<sup>2</sup> burst; thus, Schedule J was omitted

Table 18-4: ASTM D4169-16 Cycle 13 (Less Schedule I and J)

| Sequence | Temperature          | Specific Test Condition                                                          | Schedule |
|----------|----------------------|----------------------------------------------------------------------------------|----------|
| 1        | Manual Handling      | N/A                                                                              | A        |
| 2        | Vehicle Stacking     | Maximum stack height in transit vehicle                                          | C        |
| 3        | Loose Load Vibration | N/A                                                                              | F        |
| 4        | Vehicle Vibration    | Expose the packages to a 60-minute truck test followed by a 120-minute air test. | E        |
| 5        | Manual Handling      | N/A                                                                              | A        |

Anesthetic Conditioning

This conditioning was only done for select tests. Samples were conditioned for six hours with vaporized Isoflurane at 1.3% concentration and vaporized Sevoflurane at 2.6% concentration. For additional detail on the procedure used to condition samples

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with anesthetics, refer to the build and conditioning protocol and report.

### 3 Balance Coverage

#### 3.1 Testing Summary

The intent of the coverage test is to capture any patchiness of coating, at a macro-scale, on the coated surface such that it verifies coating uniformity. Coating uniformity, in turn, preserves the benefits of coating.

The principle of the coverage test is that the cationic dye (Toluidine Blue) binds to anionic substrate (due to the presence of the SO<sub>3</sub> group in Balance polymer B) and goes through metachromatic reaction and the original dye color of blue shifts to violet/purple when compared to an uncoated surface. The uncoated surface or patchiness of coating is clearly visible to an observer. This is a comparative test, so an uncoated surface must be used as a control. Toluidine Blue dye solution was prepared. Then the test and control devices were stained with the toluidine dye solution. The coverage was visually assessed by comparing the color of the test device to that of a control device. All test devices were compared to at least one uncoated control device for the coverage assessment. If the coverage was acceptable, then a "pass" was noted in the data collection sheet.

#### 3.2 Results

The t = 0 (unaged) and t = 2 YR AA sample populations of the Affinity NT Oxygenator with Balance Biosurface built with Spiderman hollow microporous polypropylene fiber produced by Medtronic Danvers satisfied the coverage design verification performance requirements.

Table 18-5 provides a summary of the results and conclusions.

Table 18-5: Balance Coverage Testing

| Requirement      | Acceptance Criteria                                                                                 | %C / %R (Data Type) | Population  | Device     | n  | Result                                   | PASS/FAIL |
|------------------|-----------------------------------------------------------------------------------------------------|---------------------|-------------|------------|----|------------------------------------------|-----------|
| Coverage (D1001) | Devices shall exhibit a metachromatic shift from blue to violet across the entire inspectable area. | 95/90 (attribute)   | t = 0       | Balance NT | 29 | 29 of 29 exhibited a metachromatic shift | PASS      |
|                  |                                                                                                     |                     | t = 2 YR AA | Balance NT | 29 | 29 of 29 exhibited a metachromatic shift | PASS      |

Abbreviations: %C = Percent Confidence, %R = Percent Reliability

### 4 Leaching

#### 4.1 Testing Summary

To verify that the Balance coating is not leaching during device usage, leaching of fluorescently labeled Polymer B, which is the top most layer of the coating, is detected and quantified. If the level of detected fluorescent signal that quantifies the potentially leached

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polymer below the specified performance requirement level, then the coating is considered non-leaching.

Theoretically, verification of Polymer B leaching also covers for Polymer A leaching because if Polymer A and B leach as bonded together, the fluorescent labeling on Polymer B shall still be detectable. If Polymer B leaches and then Polymer A leaches, the detection of Polymer B shall still be equal to or higher than Polymer A leached. Hence, leaching of both polymers are verified.

- The leachable polymer extraction circuit was set up.
- Circuit was primed with water with appropriate extraction volume required

- Water was warmed to 37°C.

Water (extraction medium) was circulated through the device for 6 hours at the maximum rated flow rate of 7 LPM while keeping the water at 37°C.

- Extract sample was collected at the end of the 6 hours.

- Fluorescent signal was read using a spectrophotometer and the leached polymer concentration was quantified.

##### 4.2 Results

The  $t = 0$  (unaged) and  $t = 2$  YR AA sample populations of the Affinity NT Oxygenator with Balance Biosurface built with Spiderman hollow microporous polypropylene fiber produced by Medtronic Danvers satisfied the Balance leaching design verification performance requirements. Table 18-6 provides a summary of the results and conclusions.

Table 18-6: Balance Leaching Testing

| Test Description | Acceptance Criteria                                                                  | %C / %R (Data Type) | Population    | Device     | n  | Result      | PASS/FAIL |
|------------------|--------------------------------------------------------------------------------------|---------------------|---------------|------------|----|-------------|-----------|
| Leaching (DI002) | Devices shall not leach fluorescently tagged polymer B at a level $\geq 0.008\%$ w/v | 95/90 (variable)    | Unaged        | Balance NT | 11 | UTL = .000% | PASS      |
|                  |                                                                                      |                     | $t = 2$ YR AA | Balance NT | 9  | UTL = .001% | PASS      |

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##### 4.3 Deviations

T=2yrAA sample size

Only 12 accelerated aged devices were coated with the fluorescent Polymer B for the leaching test. Three of these 12 devices had compromised mandrel joints and were removed from the test.

Therefore, only 9 viable test samples were available from the T=2yrAA leaching build population, which is less than the protocol requirement of  $n = 11$ . The minimum calculated sample size specified by protocol was  $n = 2$ . The minimum sample size recalculated using the T=2yrAA data was also calculated to be  $n =$

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2. The sample size used for testing exceeds both minimum sample size calculations; thus, this deviation had no impact on the outcome of the test.

## 5 Initial Gas Transfer

### Testing Summary

#### Initial Gas Transfer

The testing procedure is consistent with ISO and FDA guidelines and describes blood preparation, test apparatus, device preparation, data gathering, and data analysis methods for oxygen ( $O_2$ ) and carbon dioxide ( $CO_2$ ) transfer. Per the procedure, Bovine blood was conditioned to venous conditions and adjusted to  $12 \pm 0.3$  tHb. The gas outlet vents of all test and control devices were sealed (to ensure gas was expelled only from the gas outlet port where they are measured).

Devices were primed and de-aired, and then placed into the test circuit one at a time. Blood and gas readings were recorded for each specified time point and were used to calculate the Oxygen ( $O_2$ ) and Carbon Dioxide ( $CO_2$ ) transfer rates. In addition, blood and gas inlet and outlet pressures were recorded to determine the blood and gas path pressure drops.

#### Test Results and Conclusion

The  $t = 0$  (unaged) and  $t = 2$  YR AA sample populations of the Affinity NT Oxygenator with Balance Biosurface built with Spiderman hollow microporous polypropylene fiber produced by Medtronic Danvers satisfied the initial gas transfer design verification performance requirements. Table 18-7 provides a summary of the results and conclusions.

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| Table 10-7: Initial Gas Transfer Testing                                                                                                                                                                                                                                                                                                                |                                                 |                |               |    |                     |                   |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------|----------------|---------------|----|---------------------|-------------------|
| Requirements &<br>Acceptance Criteria                                                                                                                                                                                                                                                                                                                   | %C/%H<br>(data type)                            | Population     | Device        | n  | Result              | PASS/FAIL         |
| D0003<br>O <sub>2</sub> Transfer LTL ≥ 354<br>mL/min                                                                                                                                                                                                                                                                                                    | 95/90<br>(variable)                             | t = 0          | Balance<br>NT | 11 | LTL = 377<br>mL/min | PASS              |
|                                                                                                                                                                                                                                                                                                                                                         |                                                 | t = 2 YR<br>AA | Balance<br>NT | 11 | LTL = 375<br>mL/min | PASS              |
| D0004<br>P value is NOT less than<br>0.05 for a one-sided 2<br>sample t-test comparing<br>the gas transfer rates of<br>the test and predicate<br>samples with utilization<br>of an allowable margin.<br>If the test device mean is<br>lower than the predicate<br>device mean<br>OR<br>Test device mean is<br>greater than the predicate<br>device mean | 95/90 <sup>1</sup><br>(variable<br>comparative) | t = 0          | Balance<br>NT | 11 | 395 ± 8 mL/min      | p = 0,999<br>PASS |
|                                                                                                                                                                                                                                                                                                                                                         |                                                 |                | Tritium<br>NT | 11 | 400 ± 6 mL/min      |                   |
|                                                                                                                                                                                                                                                                                                                                                         |                                                 | t = 2 YR<br>AA | Balance<br>NT | 11 | 395 ± 8 mL/min      | p = 0,999<br>PASS |
|                                                                                                                                                                                                                                                                                                                                                         |                                                 |                | Tritium<br>NT | 11 | 396 ± 6 mL/min      |                   |
|                                                                                                                                                                                                                                                                                                                                                         |                                                 | t = 0          | Balance<br>NT | 11 | LTL = 279<br>mL/min | PASS <sup>2</sup> |
|                                                                                                                                                                                                                                                                                                                                                         |                                                 |                | Balance<br>NT | 11 | LTL = 274<br>mL/min |                   |
| D0005<br>CO <sub>2</sub> Transfer LTL ≥ 204<br>mL/min                                                                                                                                                                                                                                                                                                   | 95/90 <sup>1</sup><br>(variable<br>comparative) | t = 0          | Balance<br>NT | 11 | 200 ± 8 mL/min      | p = 0,999<br>PASS |
|                                                                                                                                                                                                                                                                                                                                                         |                                                 |                | Tritium<br>NT | 11 | 302 ± 5 mL/min      |                   |
|                                                                                                                                                                                                                                                                                                                                                         |                                                 | t = 2 YR<br>AA | Balance<br>NT | 11 | 200 ± 8 mL/min      | p = 0,999<br>PASS |
|                                                                                                                                                                                                                                                                                                                                                         |                                                 |                | Tritium<br>NT | 11 | 300 ± 6 mL/min      |                   |
|                                                                                                                                                                                                                                                                                                                                                         |                                                 | t = 0          | Balance<br>NT | 11 | UTL = 127<br>mmHg   | PASS              |
|                                                                                                                                                                                                                                                                                                                                                         |                                                 |                | Balance<br>NT | 11 | UTL = 123<br>mmHg   |                   |

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| Requirement ID & Acceptance Criteria                                                                                                                                                                                    | %C/%R (data type)   | Population     | Device        | n  | Result                                                                                                                                                                 | PASS/FAIL         |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|---------------|----|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------|
| D1008<br>Gas side pressure drop UTL at 2.0 LPM gas flow ≤ blood side pressure drop LTL at 1.0 LPM blood flow; and, Gas side pressure drop UTL at 14.0 LPM gas flow ≤ blood side pressure drop LTL at 7.0 LPM blood flow | 95/90<br>(variable) | t = 0          | Balance<br>NT | 11 | Min Flow<br>UTL <sub>gas</sub> = 3<br>mmHg<br>LTL <sub>blood</sub> = 7<br>mmHg<br><br>Max Flow<br>UTL <sub>gas</sub> = 23<br>mmHg<br>LTL <sub>blood</sub> = 95<br>mmHg | PASS              |
|                                                                                                                                                                                                                         |                     | t = 2 YR<br>AA | Balance<br>NT | 11 | Min Flow<br>UTL <sub>gas</sub> = 4<br>mmHg<br>LTL <sub>blood</sub> = 8<br>mmHg<br><br>Max Flow<br>UTL <sub>gas</sub> = 24<br>mmHg<br>LTL <sub>blood</sub> = 89<br>mmHg | PASS              |
| D1037<br>Sampling Port Accuracy LTL ≥ -10% and UTL ≤ +10%                                                                                                                                                               | 95/90<br>(variable) | t = 0          | Balance<br>NT | 11 | LTL = -3.0%<br>UTL = 4.9%                                                                                                                                              | PASS <sup>1</sup> |
|                                                                                                                                                                                                                         |                     | t = 2 YR<br>AA | Balance<br>NT | 11 | LTL = -7.1%<br>UTL = 8.0%                                                                                                                                              | PASS              |

Abbreviations: %C = percent confidence, %R = percent reliability, LTL = lower tolerance limit, UTL = upper tolerance limit  
<sup>1</sup> Due to the comparative nature of this test, a reliability value cannot be assigned.  
<sup>2</sup> See section on Test Failures below.

## 6 Time Dependent Gas Transfer

### 6.1 Testing Summary

Unlike the Initial Gas Transfer test, the intent of this time dependent (six hour) gas transfer test is to measure the oxygenator's gas exchange over a six-hour time period.

The testing procedure is consistent with ISO and FDA guidelines and describes blood preparation, test apparatus, device preparation, data gathering, and data analysis methods for oxygen (O<sub>2</sub>) and carbon dioxide (CO<sub>2</sub>) transfer tests. Per the procedure, Bovine blood was conditioned to venous conditions and adjusted to 12 +/- 0.3 tHb. The gas outlet vents of all test and control devices were sealed (to force all expelled gases out of only the gas outlet port where they can be measured).

In this test, oxygen and carbon dioxide transfer was evaluated at multiple time points over a sixhour time period. In between these blood draw intervals, devices were connected to a recirculation loop.

The test was run with four circuits set up at a time. Each circuit was connected to its own recirculation loop intended for recirculation of blood between gas transfer evaluation time points. All four circuits shared the same venous blood reservoir and were connected via the gas transfer loop. Blood and gas readings were recorded for each specified time point and were used to calculate the Oxygen (O<sub>2</sub>) and Carbon Dioxide (CO<sub>2</sub>) transfer rates. In addition, blood and gas inlet and outlet pressures were recorded to determine the blood and gas path pressure drops.

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## 6.2 Test Results and Conclusion

The t = 0 (unaged) and t = 2 YR AA sample populations of the Affinity NT Oxygenator with Balance Biosurface built with Spiderman hollow microporous polypropylene fiber produced by Medtronic Danvers satisfied the time dependent gas transfer design verification performance requirements. Table 18-8 provides a summary of the results and conclusions.

Table 18-8: Time Dependent Gas Transfer Testing

| Requirement ID & Acceptance Criteria                                                                                                                                                                                                                                                                                                                           | %C / %R (Data Type)                          | Population  | Device                  | n | Result ( $\mu \pm \sigma$ ) | PASS/FAIL        |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------|-------------|-------------------------|---|-----------------------------|------------------|
| <b>Oxygen Transfer (D1004)</b><br>P value is NOT less than 0.05 for a one-sided 2 samples t-test comparing the six-hour gas transfer rates of the test and predicate samples with utilization of an allowable margin, if the test device mean is lower than the predicate device mean OR<br>Test device mean is greater than the predicate device mean         | 95/NA <sup>1</sup><br>(variable comparative) | t = 0       | Balance NT              | 6 | .990 $\pm$ .022             | p = .588<br>PASS |
|                                                                                                                                                                                                                                                                                                                                                                |                                              |             | Trillium NT (predicate) | 6 | 1.000 $\pm$ .030            |                  |
|                                                                                                                                                                                                                                                                                                                                                                |                                              | t = 2 YR AA | Balance NT              | 6 | .994 $\pm$ .016             | p = .790<br>PASS |
|                                                                                                                                                                                                                                                                                                                                                                |                                              |             | Trillium NT (predicate) | 6 | 1.000 $\pm$ .015            |                  |
|                                                                                                                                                                                                                                                                                                                                                                |                                              | t = 0       | Balance NT              | 6 | .966 $\pm$ .011             | p = .100<br>PASS |
|                                                                                                                                                                                                                                                                                                                                                                |                                              |             | Trillium NT (Predicate) | 6 | 1.000 $\pm$ .017            |                  |
| <b>Carbon Dioxide Transfer (D1005)</b><br>P value is NOT less than 0.05 for a one-sided 2 samples t-test comparing the six-hour gas transfer rates of the test and predicate samples with utilization of an allowable margin, if the test device mean is lower than the predicate device mean OR<br>Test device mean is greater than the predicate device mean | 95/NA <sup>1</sup><br>(variable comparative) | t = 2 YR AA | Balance NT              | 6 | .994 $\pm$ .027             | p = .779<br>PASS |
|                                                                                                                                                                                                                                                                                                                                                                |                                              |             | Trillium NT (Predicate) | 6 | 1.000 $\pm$ .044            |                  |
|                                                                                                                                                                                                                                                                                                                                                                |                                              |             |                         |   |                             |                  |

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| Requirement ID & Acceptance Criteria                                                                                                                                                                                         | %C / %R (Data Type)                      | Population  | Device     | n | Result (μ ± σ)                           | PASS/FAIL                                |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------|-------------|------------|---|------------------------------------------|------------------------------------------|
| FDA Gas Transfer (FDA Guidance)<br>Collect gas transfer data at 90-minute increments over six hours at multiple blood flow rates (minimum, nominal, maximum) and blood:gas ratios (1:1, 1:2, 1:0.5) at 100% FiO <sub>2</sub> | N/A - characterization only <sup>1</sup> | t = 0       | Balance NT | 6 | See Report 10845429DOC (Attachment 18-8) | N/A - characterization only <sup>1</sup> |
|                                                                                                                                                                                                                              |                                          | t = 2 YR AA | Balance NT | 6 | See Report 10845429DOC (Attachment 18-8) | N/A - characterization only <sup>1</sup> |

Abbreviations: %C = percent confidence; %R = percent reliability

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<sup>2</sup> This test is not evaluated against a specification. Testing was conducted to provide FDA with data per the FDA guidance

## 7 Prime Volume

### 7.1 Background

#### Prime Volume

A low prime volume in the cardiopulmonary bypass circuit is desirable. As the perfusion circuit prime volume decreases, the patient hemodilution also decreases. The device dry weight was measured. The device was filled with test media and the wet weight measured. The prime volume was calculated using the following equation:

$$\text{Prime Volume} = \frac{(\text{Wet Weight} - \text{Dry Weight})}{1.0046}$$

#### Residual Prime Volume

The residual prime volume is measured using saline but can be thought of as the amount of blood lost from the patient – this is the volume of fluid that remains in the device following the bypass procedure. A low residual prime volume is also desirable.

Following prime volume testing, the device was held such that the blood inlet and outlet were facing downwards. The device was drained for 20 seconds and then clamped off. The device was weighed again, and the device residual prime weight was recorded. The residual prime volume was calculated using the following equation:

$$\text{Residual Prime Volume} = \frac{(\text{Residual Prime Weight} - \text{Dry Prime Weight})}{1.0046}$$

The above equations incorporate the conversion factor 1 mL saline = 1.0046 g.

### 7.2 Test Results and Conclusion

The t = 0 (unaged) and t = 2 YR AA sample populations of the Affinity NT Oxygenator with Balance Biosurface built with

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Spiderman hollow microporous polypropylene fiber produced by Medtronic Danvers satisfied the prime volume design verification performance requirements. Table 18-9 provides a summary of the results and conclusions.

Table 18-9: Prime Volume Testing

| Design Input Requirement (ID) | Specification Acceptance Criteria                                                                                                                                                                              | %Confidence / %Reliability (Data Type) | Population  | Device     | n  | Result (μ ± σ) | PASS/FAIL                            |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------|-------------|------------|----|----------------|--------------------------------------|
| Prime Volume (DI009)          | Prime Volume LTL ≥ 250 mL<br>Prime Volume UTL ≤ 302 mL                                                                                                                                                         | 90/90 (variable)                       | Unaged      | Balance NT | 11 | 280 ± 2 mL     | LTL = 275 mL<br>UTL = 285 mL<br>PASS |
|                               |                                                                                                                                                                                                                |                                        | t = 2 YR AA | Balance NT | 11 | 279 ± 2 mL     | LTL = 274 mL<br>UTL = 284 mL<br>PASS |
| Residual Prime Volume (DI010) | n/a – characterization only. Specifications and acceptance criteria do not exist for residual prime volume. The testing is being completed to support residual prime volume claims made within the device IFU. | n/a – characterization only            | Unaged      | Balance NT | 11 | 149 ± 9 mL     | n/a – characterization only          |
|                               |                                                                                                                                                                                                                |                                        | t = 2 YR AA | Balance NT | 11 | 129 ± 6 mL     | n/a – characterization only          |

## 8 Maximum Blood Flow Trauma Testing Summary

The test procedure is consistent with ISO and FDA guidelines and describes blood preparation, test apparatus, device preparation, data gathering, and data analysis methods for Hemolysis, White Blood Cell, Platelet and Platelet Function tests. Per the procedure, conditioned Bovine blood was utilized as the test media. Platelet works assay was done at time zero, and at 360 minutes on each device. Platelet assays were run in triplicates using a blood sample from a single draw.

The intent of the Blood Trauma test is to measure the damage of the blood over time. Therefore, the blood conditions of each test circuit were required to fall within the FDA required conditions at the first circuit draw time point as described in Table 18-7.

Table 18-7: Sample FDA Required Blood Conditions

| Blood Factor                     | Level                                      |
|----------------------------------|--------------------------------------------|
| Hemoglobin Concentration         | 12 ± 1 g/dL <sup>2</sup>                   |
| pCO <sub>2</sub>                 | 40 ± 5 mmHg                                |
| PH                               | 7.4 ± 0.1                                  |
| Activated Clotting Time          | ≥ 300 s                                    |
| Blood Glucose                    | 100-300 mg/dL                              |
| Circuit Temperature <sup>1</sup> | 37 ± 2°C (ISO 7199: 37 ± 1°C) <sup>1</sup> |

<sup>1</sup> Temperature shall be maintained over the course of the six-hour test. ISO 7199 requires circuit temperature be maintained at 37 ± 1°C. This is more conservative than the FDA requirements, thus the ISO requirements shall be maintained during the test.

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The devices were then prepared for thrombus inspection after completion of the Blood Trauma Test, using following steps:

- Blood was filtered from each circuit using a 38 μm Affinity Filter
- Device was gently removed from circuit and rinsed with 0.9% saline. The flow of saline was low and just enough to rinse the device. Extra care was taken to not disrupt results.
- All excess saline was gravity drained

Both test, predicate and reference devices (and circuits) were inspected for the presence of thrombus. The visual inspection was gross in nature and data was recorded as "yes" or "no" in terms of presence of thrombus.

#### Testing Results and Conclusion

The t = 0 (unaged) and t = 2 YR AA sample populations of the Affinity NT Oxygenator with Balance Biosurface built with Spiderman hollow microporous polypropylene fiber produced by Medtronic Danvers satisfied the maximum flow blood trauma design verification performance requirements. Table 18-8 provides a summary of the results and conclusions.

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Table 18-8: Maximum Blood Flow Trauma Testing

| Design Input Requirement (ID)            | Specification Acceptance Criteria                                                                                                                                                                                                                                                                                      | %Confidence / %Reliability (Data Type) | Population    | Device                     | n | Result ( $\bar{x} \pm s$ ) | PASS/FAIL                          |
|------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------|---------------|----------------------------|---|----------------------------|------------------------------------|
| Hemolysis (MIH) (DI011)                  | P value is NOT less than 0.05 for a one-sided two samples t-test comparing 6-hour test device MIH normalized counts and predicate device MIH, with utilization of an allowable margin, if the test device mean is greater than the predicate device mean, OR test device mean is less than the predicate device mean.  | 95/NA* (variable comparative)          | $t = 0$       | Balance NT (Test)          | 6 | 1.262 $\pm$ .519           | $p = .686$<br>PASS                 |
|                                          |                                                                                                                                                                                                                                                                                                                        |                                        |               | Trillium NT (Predicate)    | 6 | 1.000 $\pm$ .165           |                                    |
|                                          |                                                                                                                                                                                                                                                                                                                        |                                        |               | Balance Fusion (Reference) | 6 | 1.592 $\pm$ .462           |                                    |
|                                          |                                                                                                                                                                                                                                                                                                                        |                                        | $t = 2$ YR AA | Balance NT (Test)          | 6 | 1.085 $\pm$ .505           | $p = .821$<br>PASS                 |
|                                          |                                                                                                                                                                                                                                                                                                                        |                                        |               | Trillium NT (Predicate)    | 6 | 1.000 $\pm$ .537           |                                    |
|                                          |                                                                                                                                                                                                                                                                                                                        |                                        |               | Balance Fusion (Reference) | 6 | 1.505 $\pm$ .486           |                                    |
| White Blood Cell Retention (WBC) (DI012) | P value is NOT less than 0.05 for a one-sided two samples t-test comparing 6-hour test device WBC normalized counts and predicate device WBC, with utilization of an allowable margin, if the test device mean is lower than the predicate device mean, OR test device mean is greater than the predicate device mean. | 95/NA* (variable comparative)          | $t = 0$       | Balance NT (Test)          | 6 | .989 $\pm$ .080            | $p = 1.000$<br>PASS                |
|                                          |                                                                                                                                                                                                                                                                                                                        |                                        |               | Trillium NT (Predicate)    | 6 | 1.000 $\pm$ .027           |                                    |
|                                          |                                                                                                                                                                                                                                                                                                                        |                                        |               | Balance Fusion (Reference) | 6 | 1.006 $\pm$ .059           |                                    |
|                                          |                                                                                                                                                                                                                                                                                                                        |                                        | $t = 2$ YR AA | Balance NT (Test)          | 6 | 1.432 $\pm$ .632           | $p_{TEST} > p_{PREDICATE}$<br>PASS |
|                                          |                                                                                                                                                                                                                                                                                                                        |                                        |               | Trillium NT (Predicate)    | 6 | 1.000 $\pm$ .200           |                                    |
|                                          |                                                                                                                                                                                                                                                                                                                        |                                        |               | Balance Fusion (Reference) | 6 | .971 $\pm$ .130            |                                    |

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| Design Input Requirement (ID)    | Specification Acceptance Criteria                                                                                                                                                                                                                                                                                                            | %Confidence / %Reliability (Data Type) | Population  | Device                     | n | Result ( $\mu \pm \sigma$ ) | PASS/FAIL                  |
|----------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------|-------------|----------------------------|---|-----------------------------|----------------------------|
| Platelet (PLT) Retention (DI013) | P value is NOT less than 0.05 for a one-sided two samples t-test comparing 6-hour test device PLT normalized counts and predicate device PLT, with utilization of an allowable margin, if the test device mean is lower than the predicate device mean, OR Test device mean is greater than the predicate device mean.                       | 95/NA* (variable comparative)          | t = 0       | Balance NT (Test)          | 6 | 1.070 $\pm$ .085            | μTEST > μPREDICATE<br>PASS |
|                                  |                                                                                                                                                                                                                                                                                                                                              |                                        |             | Trillium NT (Predicate)    | 6 | 1.000 $\pm$ .019            |                            |
|                                  |                                                                                                                                                                                                                                                                                                                                              |                                        |             | Balance Fusion (Reference) | 6 | .998 $\pm$ .042             |                            |
|                                  |                                                                                                                                                                                                                                                                                                                                              |                                        | t = 2 YR AA | Balance NT (Test)          | 6 | 1.058 $\pm$ .055            | μTEST > μPREDICATE<br>PASS |
|                                  |                                                                                                                                                                                                                                                                                                                                              |                                        |             | Trillium NT (Predicate)    | 6 | 1.000 $\pm$ .026            |                            |
|                                  |                                                                                                                                                                                                                                                                                                                                              |                                        |             | Balance Fusion (Reference) | 6 | 1.012 $\pm$ .061            |                            |
| Function Platelets (DI014)       | P value is NOT less than 0.05 for a one-sided two samples t-test comparing 6-hour test device functional PLT normalized counts and predicate device functional PLT, with utilization of an allowable margin, if the test device mean is lower than the predicate device mean, OR test device mean is greater than the predicate device mean. | 95/NA* (variable comparative)          | t = 0       | Balance NT (Test)          | 6 | 1.087 $\pm$ .152            | μTEST > μPREDICATE<br>PASS |
|                                  |                                                                                                                                                                                                                                                                                                                                              |                                        |             | Trillium NT (Predicate)    | 6 | 1.000 $\pm$ .043            |                            |
|                                  |                                                                                                                                                                                                                                                                                                                                              |                                        |             | Balance Fusion (Reference) | 6 | .897 $\pm$ .119             |                            |
|                                  |                                                                                                                                                                                                                                                                                                                                              |                                        | t = 2 YR AA | Balance NT (Test)          | 6 | 1.223 $\pm$ .102            | μTEST > μPREDICATE<br>PASS |
|                                  |                                                                                                                                                                                                                                                                                                                                              |                                        |             | Trillium NT (Predicate)    | 6 | 1.000 $\pm$ .117            |                            |
|                                  |                                                                                                                                                                                                                                                                                                                                              |                                        |             | Balance Fusion (Reference) | 6 | 0.961 $\pm$ .162            |                            |

| Design Input Requirement (ID) | Specification Acceptance Criteria                                                                                                                                                           | %Confidence / %Reliability (Data Type) | Population  | Device                     | n | Result ( $\mu \pm \sigma$ ) | PASS/FAIL                   |
|-------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------|-------------|----------------------------|---|-----------------------------|-----------------------------|
| Thrombus Inspection (DI016)   | N/A = characterization only. Specifications and acceptance criteria do not exist for thrombus inspection. The testing is being completed to provide the FDA with data per the FDA guidance. | NA = characterization only             | Unaged      | Balance NT (Test)          | 6 | Thrombus present            | N/A = characterization only |
|                               |                                                                                                                                                                                             |                                        |             | Trillium NT (Predicate)    | 6 | Thrombus present            |                             |
|                               |                                                                                                                                                                                             |                                        |             | Balance Fusion (Reference) | 6 | Thrombus present            |                             |
|                               |                                                                                                                                                                                             |                                        | t = 2 YR AA | Balance NT (Test)          | 6 | Thrombus present            | N/A = characterization only |
|                               |                                                                                                                                                                                             |                                        |             | Trillium NT (Predicate)    | 6 | Thrombus present            |                             |
|                               |                                                                                                                                                                                             |                                        |             | Balance Fusion (Reference) | 6 | Thrombus present            |                             |

\*Due to the comparative nature of this test, a reliability value cannot be assigned

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## 9 Minimum Blood Flow Trauma

### Testing Results and Conclusion

The t = 0 (unaged) and t = 2 YR AA sample populations of the Affinity NT Oxygenator with Balance Biosurface built with Spiderman hollow microporous polypropylene fiber produced by Medtronic Danvers satisfied the minimum flow blood trauma design verification performance requirements. Table 18-9 provides a summary of the results and conclusions.

Table 18-9: Minimum Blood Flow Trauma Testing

| Design Input Requirement (ID) | Specification Acceptance Criteria                                                                                                                                                                                                                                                                                                     | %Confidence / %Reliability (Data Type)    | Population  | Device                     | n | Result (μ ± σ)   | PASS/FAIL                    |
|-------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------|-------------|----------------------------|---|------------------|------------------------------|
| Functional Platelets (D1015)  | P value is NOT less than 0.05 for a one-sided 2 samples t-test comparing the Functional Platelets of the test and predicate samples with utilization of an allowable margin, if the test device mean is less than the predicate device mean OR Test device mean is greater than the predicate device mean                             | 95/NA <sup>1</sup> (variable comparative) | t = 0       | Balance NT (Test)          | 6 | .910 ± .104      | p = .590<br>PASS             |
|                               |                                                                                                                                                                                                                                                                                                                                       |                                           |             | Trillium NT (Predicate)    | 6 | 1.000 ± .051     |                              |
|                               |                                                                                                                                                                                                                                                                                                                                       |                                           |             | Balance Fusion (Reference) | 6 | .951 ± .167      |                              |
|                               |                                                                                                                                                                                                                                                                                                                                       |                                           | t = 2 YR AA | Balance NT (Test)          | 6 | 1.038 ± .149     | P TEST > P PREDICATE<br>PASS |
|                               |                                                                                                                                                                                                                                                                                                                                       |                                           |             | Trillium NT (Predicate)    | 6 | 1.000 ± .148     |                              |
|                               |                                                                                                                                                                                                                                                                                                                                       |                                           |             | Balance Fusion (Reference) | 6 | 0.892 ± .203     |                              |
| Thrombus Inspection (D1017)   | N/A – characterization only – Devices shall be inspected for thrombus after six hours of blood circulation under the minimum rated flow and results shall be recorded. Specifications and acceptance criteria do not exist for thrombus inspection. The testing is being completed to provide the FDA with data per the FDA guidance. | NA – characterization only                | Unaged      | Balance NT (Test)          | 6 | No Thrombus      | N/A – characterization only  |
|                               |                                                                                                                                                                                                                                                                                                                                       |                                           |             | Trillium NT (Predicate)    | 6 | Thrombus present |                              |
|                               |                                                                                                                                                                                                                                                                                                                                       |                                           |             | Balance Fusion (Reference) | 6 | Thrombus present |                              |
|                               |                                                                                                                                                                                                                                                                                                                                       |                                           | t = 2 YR AA | Balance NT (Test)          | 6 | Thrombus present | N/A – characterization only  |
|                               |                                                                                                                                                                                                                                                                                                                                       |                                           |             | Trillium NT (Predicate)    | 6 | Thrombus present |                              |
|                               |                                                                                                                                                                                                                                                                                                                                       |                                           |             | Balance Fusion (Reference) | 6 | Thrombus present |                              |

<sup>1</sup> Due to the comparative nature of this test, a reliability value cannot be assigned

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The intent of this test is to test the device port strength. The device port strength must be greater than the force the ports may experience during foreseeable use.

A port break arm was attached to the Instron load cell, and the oxygenator under test was fixed onto the Instron such that the port axis is normal to the platform and concentric with the end of the port break arm. The port was compressed at a rate of 6 inches per minute. The maximum load was recorded.

## 10.2 Testing Results and Conclusion

The t = 0 (unaged) and t = 2 YR AA sample populations of the Affinity NT Oxygenator with Balance Biosurface built with Spiderman hollow microporous polypropylene fiber produced by Medtronic Danvers satisfied the port break design verification performance requirements. Table 18-10 provides a summary of the results and conclusions.

Table 18-10: Port Break Testing

| Design Input Requirement (ID)    | Specification Acceptance Criteria | % Confidence / % Reliability (Data Type) | Population  | Port               | n  | Result (μ ± σ) | PASS/FAIL             |
|----------------------------------|-----------------------------------|------------------------------------------|-------------|--------------------|----|----------------|-----------------------|
| Blood Port Break (D1026)         | Blood Port Break LTL ≥ 22.6 lbf   | 95/90 (variable)                         | Unaged      | Blood Inlet        | 11 | 226 ± 8 lbf    | LTL = 207 lbf<br>PASS |
|                                  |                                   |                                          |             | Blood Outlet       | 11 | 161 ± 7 lbf    | LTL = 146 lbf<br>PASS |
|                                  |                                   |                                          | t = 2 YR AA | Blood Inlet        | 11 | 210 ± 5 lbf    | LTL = 199 lbf<br>PASS |
|                                  |                                   |                                          |             | Blood Outlet       | 11 | 152 ± 2 lbf    | LTL = 147 lbf<br>PASS |
| Water Port Break (D1027)         | Water Port Break LTL ≥ 22.6 lbf   | 95/90 (variable)                         | Unaged      | Water Inlet        | 11 | 360 ± 3 lbf    | LTL = 353 lbf<br>PASS |
|                                  |                                   |                                          |             | Water Outlet       | 11 | 355 ± 4 lbf    | LTL = 345 lbf<br>PASS |
|                                  |                                   |                                          | t = 2 YR AA | Water Inlet        | 11 | 353 ± 2 lbf    | LTL = 348 lbf<br>PASS |
|                                  |                                   |                                          |             | Water Outlet       | 11 | 342 ± 11 lbf   | LTL = 318 lbf<br>PASS |
| Gas Inlet Port Break (D1028)     | Gas Port Break LTL ≥ 22.6 lbf     | 95/90 (variable)                         | Unaged      | Gas Inlet          | 11 | 161 ± 6 lbf    | LTL = 149 lbf<br>PASS |
|                                  |                                   |                                          | t = 2 YR AA | Gas Inlet          | 11 | 130 ± 14 lbf   | LTL = 99 lbf<br>PASS  |
| Recirculation Port Break (D1029) | Recirc Port Break LTL ≥ 22.6 lbf  | 95/90 (variable)                         | Unaged      | Recirculation Port | 11 | 125 ± 3 lbf    | LTL = 119 lbf<br>PASS |
|                                  |                                   |                                          | t = 2 YR AA | Recirculation Port | 11 | 121 ± 4 lbf    | LTL = 112 lbf<br>PASS |

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The at to easily place tubing facilitates setup of the CPB circuit. The tubing must also remain attached to ports under forces experienced during foreseeable use.

The port under test was prepared as per the test method. The tubing was first pulled at a constant force of 15 N (3.4 lbf) for 15 seconds. After 15 seconds, the tubing was pulled at 5 in/min until failure. The maximum load at failure was recorded by Instron program.

## 11.2 Testing Results and Conclusion

The t = 0 (unaged) and t = 2 YR AA sample populations of the Affinity NT Oxygenator with Balance Biosurface built with Spiderman hollow microporous polypropylene fiber produced by Medtronic Danvers satisfied the tubing pull design verification performance requirements.

Table 18-11 provides a summary of the results and conclusions.

Table 18-11: Tubing Pull Testing

| Design Input Requirement (ID)       | Acceptance Criteria                                      | %Confidence / %Reliability (Data Type) | Population  | Port          | n  | Result (μ ± σ)       | PASS/FAIL             |
|-------------------------------------|----------------------------------------------------------|----------------------------------------|-------------|---------------|----|----------------------|-----------------------|
| Tubing Pull (Removal Force) (DI031) | Tubing Pull Off Force LTL ≥ 22 lbf                       | 95/90 (variable)                       | Unaged      | Blood Inlet   | 11 | 159 ± 7 lbf          | LTL = 142 lbf<br>PASS |
|                                     |                                                          |                                        |             | Blood Outlet  | 11 | 150 ± 5 lbf          | LTL = 138 lbf<br>PASS |
|                                     |                                                          |                                        |             | Recirculation | 11 | 105 ± 3 lbf          | LTL = 99 lbf<br>PASS  |
|                                     |                                                          |                                        |             | Gas Inlet     | 11 | 107 ± 3 lbf          | LTL = 101 lbf<br>PASS |
|                                     |                                                          |                                        | t = 2 YR AA | Blood Inlet   | 11 | 126 ± 8 lbf          | LTL = 109 lbf<br>PASS |
|                                     |                                                          |                                        |             | Blood Outlet  | 11 | 121 ± 8 lbf          | LTL = 103 lbf<br>PASS |
|                                     |                                                          |                                        |             | Recirculation | 11 | 83 ± 5 lbf           | LTL = 72 lbf<br>PASS  |
|                                     |                                                          |                                        |             | Gas Inlet     | 11 | 83 ± 4 lbf           | LTL = 74 lbf<br>PASS  |
| Tubing Pull (15s hold) (DI032)      | Tubing pull withstands a force of 3.4 lbf for 15 seconds | n/a - characterization only            | Unaged      | Blood Inlet   | 29 | 29 of 29 no pull off | PASS                  |
|                                     |                                                          |                                        |             | Blood Outlet  | 29 | 29 of 29 no pull off | PASS                  |
|                                     |                                                          |                                        |             | Recirculation | 29 | 29 of 29 no pull off | PASS                  |
|                                     |                                                          |                                        |             | Gas Inlet     | 22 | 22 of 22 no pull off | PASS                  |
|                                     |                                                          |                                        | t = 2 YR AA | Blood Inlet   | 29 | 29 of 29 no pull off | PASS                  |
|                                     |                                                          |                                        |             | Blood Outlet  | 29 | 29 of 29 no pull off | PASS                  |
|                                     |                                                          |                                        |             | Recirculation | 29 | 29 of 29 no pull off | PASS                  |
|                                     |                                                          |                                        |             | Gas Inlet     | 22 | 22 of 22 no pull off | PASS                  |

## 12 Temperature Measurement Adaptor (TMA) Accuracy

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## 12.1 Testing Summary

During cardiopulmonary bypass, blood is cooled to lower a patient's core temperature to reduce metabolic load. The patient must also be brought back to normal core temperature before coming off bypass. The TMA is the interface between the patient's blood and the temperature probe.

Blood temperature is monitored at the TMA for an indication of the arterial blood's temperature returning to the patient.

Prior to commencement of the test, the temperature probe and T-type thermocouple used for the test were checked to verify they are reading temperatures within 0.2°C of each other. Using saline as the test media, devices were placed into the test circuit one at a time. The TMA and bulk flow temperatures were recorded to determine the accuracy of the temperature measurement at the TMA interface.

The TMA accuracy was calculated using the temperature values recorded at the TMA and at the bulk flow thermocouple locations. The TMA accuracy is defined as:

TMA Accuracy = (TMA temperature) – (Bulk flow thermocouple temperature)

## 12.2 Testing Results and Conclusion

The Affinity NT Oxygenator with Balance Biosurface built with Spiderman hollow microporous polypropylene fiber produced by Medtronic Danvers satisfied the TMA Accuracy design verification performance requirements. Table 18-12 provides a summary of the results and conclusions.

Table 18-12: TMA Accuracy Testing

| Design Input Requirement (ID) | Acceptance Criteria                                                                      | %Confidence / %Reliability (Data Type) | Population  | Device     | n  | Result                                                   | PASS/FAIL |
|-------------------------------|------------------------------------------------------------------------------------------|----------------------------------------|-------------|------------|----|----------------------------------------------------------|-----------|
| TMA Accuracy (D1033)          | TMA Accuracy LTL $\geq -1.5^\circ\text{C}$<br>TMA Accuracy UTL $\leq +1.5^\circ\text{C}$ | 90/90 (variable)                       | Unaged      | Balance NT | 11 | LTL = $-1.4^\circ\text{C}$<br>UTL = $-0.7^\circ\text{C}$ | NA*       |
|                               |                                                                                          |                                        | t = 2 YR AA | Balance NT | 11 | LTL = $-1.0^\circ\text{C}$<br>UTL = $-0.3^\circ\text{C}$ | PASS      |

### \*Deviations

The TMA Accuracy test method is not a released standalone test method, instead the test procedure was outlined in the test protocol. Consequently, there is no system record that the test operator was trained to this embedded test method prior to 30 July 2018 (unaged) test execution and therefore, pass/fail determination was left as 'NA' for t=0. The operator was subsequently trained on the test method prior to the t = 2 YR AA test execution. TMA accuracy

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performance is governed by part geometry and is not impacted by aging. Thus, the T=2yrAA results documented in this report sufficiently demonstrates product conformance to requirements.

### 13 Dust Cap Removal Force and Particulate Inspection

#### 13.1 Testing Summary

The force required to remove the dust cap must be less than or equal to the force which may be exerted on the caps upon removal during clinical use. This test is for characterization/information purposes only because the removal force is a customer requirement and has no impact on patient safety. Immediately after the removal of the blood inlet and outlet dust caps, the blood inlet and outlet ports undergo particulate inspection.

The device to be tested was placed into a test fixture with the test port facing up. An Instron with pneumatic grips was secured to the top of the dust cap and pulled at a constant speed of 5 in/min. The maximum force required to remove the dust cap was recorded.

After dust cap removal, the surfaces of the port were inspected, and particulate observations were recorded.

#### 13.2 Testing Results and Conclusion

The t = 0 (unaged) sample population of the Affinity NT Oxygenator with Balance Biosurface built with Spiderman hollow microporous polypropylene fiber produced by Medtronic Danvers satisfied the dust cap particulate removal design verification performance requirements. Table 18-13 provides a summary of the results and conclusions.

Table 18-13: Dust Cap Removal Force and Particulate Inspection Testing

| Requirement (ID)               | Acceptance Criteria         | %C/%R (Data Type)           | Population | n  | Result (μ ± σ)                                                                                                                                                                | PASS/FAIL                   |
|--------------------------------|-----------------------------|-----------------------------|------------|----|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------|
| Dust Cap Removal Force (D1035) | n/a - characterization only | N/A - characterization Only | t = 0      | 11 | Blood Inlet: 2.5 ± 1.1 lbf<br>Blood Outlet: 3.4 ± 1.2 lbf<br>Gas Inlet: 6.5 ± 3.1 lbf<br>Recirc: 4.9 ± 3.4 lbf<br>Water Inlet: 3.4 ± 0.4 lbf<br>Water Outlet: 3.2 ± 0.8 lbf   | N/A - characterization only |
|                                |                             |                             |            | 11 | Blood Inlet: 2.4 ± 0.8 lbf<br>Blood Outlet: 3.0 ± 2.3 lbf<br>Gas Inlet: 11.7 ± 5.4 lbf<br>Recirc: 12.3 ± 2.2 lbf<br>Water Inlet: 4.5 ± 1.2 lbf<br>Water Outlet: 4.6 ± 1.0 lbf | N/A - characterization only |

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| Requirement (ID)                        | Acceptance Criteria                                                             | %C/%R (Data Type) | Population  | n  | Result (μ ± σ)          | PASS/FAIL         |
|-----------------------------------------|---------------------------------------------------------------------------------|-------------------|-------------|----|-------------------------|-------------------|
| Dust Cap Particulate Inspection (D1036) | No particulate shall be present on the device following removal of the dust cap | 90/90 (attribute) | Unaged      | 22 | 22 of 22 no particulate | PASS <sup>1</sup> |
|                                         |                                                                                 |                   | t = 2 YR AA | 22 | 22 of 22 no particulate | PASS              |

Abbreviations: %C = percent confidence, %R = percent reliability; N/A = not applicable  
<sup>1</sup> See deviations and test failures below

## 14 Particulate

### 14.1 Testing Summary

Particulate matter consists of particles that will not dissolve in solution (other than gas bubbles) and that are unintentionally present on a device or in a solution. Particulate matter can come from many sources during manufacturing. The testing is intended to verify that these particulate levels do not exceed the limits set forth in the test method.

The fluid pathway of each test article was flushed with 250 mL of particle-free water. The flush solution was collected in a clean container and vacuum filtered through a 0.8 μm membrane filter. A blank sample was also prepared using the same procedure, but without the test article.

The filters were microscopically examined and counted per USP<788> with particle size expressed as ≥10μm and ≥25 μm.

### 14.2 Testing Results and Conclusion

The Affinity NT Oxygenator with Balance Biosurface built with Spiderman hollow microporous polypropylene fiber produced by Medtronic Danvers satisfied the particulate design verification performance requirements. Table 18-14 provides a summary of the results and conclusions.

Table 18-14: Particulate Testing

| Requirement ID             | Acceptance Criteria      | %C/%R (Data Type) | Population | Device     | n | Result    | PASS/FAIL |
|----------------------------|--------------------------|-------------------|------------|------------|---|-----------|-----------|
| Particulate ≥10 μm (D1046) | Particulate UTL ≤ 13,750 | 95/90 (variable)  | Unaged     | Balance NT | 5 | UTL = 118 | PASS      |
| Particulate ≥25 μm (D1046) | Particulate UTL ≤ 1,375  | 95/90 (variable)  | Unaged     | Balance NT | 5 | UTL = 61  | PASS      |

Abbreviations: %C = percent confidence, %R = percent reliability

## 15 Water Side Burst

### 15.1 Testing Summary

Device is pressurized until failure to ensure the device is capable of withstanding 2x or more the rated pressure.

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Device was connected to a Rothenberger hydrostatic test pump, an Ashcroft digital pressure gauge, and burst chamber. The test pump was filled with water and was used to pressurize the water path of the device until failure. Due to limitations of the pump and tubing clamps, if 400 psi was reached without a failure the test was stopped and the maximum pressure reading was recorded.

### 15.2 Testing Results and Conclusion

The  $t=0$  (unaged) and  $t=2$  YR AA sample populations of the Affinity NT Oxygenator with Balance Biosurface built with Spiderman hollow microporous polypropylene fiber produced by Medtronic Danvers satisfied the water side burst design verification performance requirements. Table 18-15 provides a summary of the results and conclusions.

Table 18-15: Water Side Burst Testing

| Design Input Requirement (ID) | Acceptance Criteria                           | %Confidence / %Reliability (Data Type) | Population  | Device     | n  | Result        | PASS/FA |
|-------------------------------|-----------------------------------------------|----------------------------------------|-------------|------------|----|---------------|---------|
| Water Side Burst (DI023)      | The device water path burst LTL $\geq 60$ psi | 90/90 (variable)                       | Unaged      | Balance NT | 11 | LTL = 400 psi | PASS    |
|                               |                                               |                                        | $t=2$ YR AA | Balance NT | 11 | LTL = 394 psi | PASS    |

### 16 Blood and Gas Side Burst

#### 16.1 Testing Summary

Device is pressurized until failure to ensure the device is capable of withstanding 2x or more the rated pressure. Test devices were connected to a Rothenberger hydrostatic test pump, an Ashcroft digital pressure gauge, and burst chamber. The test pump was filled with water and was used to pressurize the desired flow path of the device until failure.

Upon burst, the maximum pressure reached as well as the failure mode were recorded.

#### 16.2 Testing Results and Conclusion

The  $t=0$  (unaged) and  $t=2$  YR AA test populations of the Affinity NT Oxygenator with Balance Biosurface built with Spiderman hollow microporous polypropylene fiber produced by Medtronic Danvers satisfied the blood and gas side burst performance requirements after 6 hours of gas path exposure to 1.3% Isoflurane and 2.6% Sevoflurane. Table 18-19 provides a summary of the results and conclusions.

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Table 18-19: Blood and Gas Side Burst Testing

| Requirement (ID)         | Acceptance Criteria        | %C / %R (Data Type) | Population  | n  | Result          | PASS/FAIL |
|--------------------------|----------------------------|---------------------|-------------|----|-----------------|-----------|
| Blood Path Burst (D1022) | Burst LTL $\geq 1500$ mmHg | 95/90 (variable)    | t = 0       | 11 | LTL = 6427 mmHg | PASS      |
|                          |                            |                     | t = 2 YR AA | 11 | LTL = 7229 mmHg | PASS      |
| Gas Path Burst (D1024)   | Burst LTL $\geq 98$ mmHg   | 90/90 (variable)    | t = 0       | 11 | LTL = 609 mmHg  | PASS      |
|                          |                            |                     | t = 2 YR AA | 11 | LTL = 740 mmHg  | PASS      |

Abbreviations: %C = percent confidence, %R = percent reliability.

<sup>1</sup> The initial t=2yrAA test population did not meet blood side burst requirements due to the presence of an isolated outlier (SN 5116096025) with blood side burst pressure at 252 mmHg (5 psi). Refer to Section 18.14 for additional details.

#### 16.3 Test Failures

##### T=2yrAA Blood Side Burst

Prior to burst testing, this outlier device (SN 5116096025) held 1125 mmHg (22 psi) for a duration of six hours during the pressure integrity test, thus demonstrating the initial integrity of the device. Further investigation pointed to a single impact type event leading to substrate damage proximal to the joint where the burst occurred. A series of directional cracks were visible over a length of 0.199", as well as deformation of the HEX housing near the side seam joint. A visual inspection was also performed on the remaining 28 devices from the pressure integrity test population, yielding 0 occurrences of HEX damage on other test samples. Therefore, this was noted as an isolated event unrelated to the design function or mechanical integrity of the joint adhesive. Medtronic believes the suspected root cause of this failure was due to drop/mishandling between tests. The type, size, and location of the damage as well as the verified initial integrity of this unit supports this conclusion. Though there is reason to omit this outlier based on the deduction that the unit changed state between tests and the number of remaining passing units sufficiently met the minimum risk based sample size, documentation of a mishandling event does not exist. Therefore, a population of another 11 devices were retested to increase statistical confidence and meet original protocol intent. All 11 samples passed upon retest.

#### 17 Water Side Integrity

##### 17.1 Testing Summary

The oxygenator is under constant fluid pressure and is responsible for heating and cooling the blood and transferring oxygen and

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carbon dioxide. There is pressure exerted on the water-side of the device from the water pump. The water-side must not leak for the oxygenator to function properly. As a result, a maximum rated operating pressure is specified for the water-side. The integrity of the water path is verified by subjecting the water path to a hydraulic pressure more than the specified maximum recommended operating pressure for 6 hours, and then observing for leakage.

The test circuit was constructed. Once the pump was primed and the circuit was debubbled, the circuit was pressurized to  $45 \pm 3$  psi and checked for leaks. The circuit pressure was held for six hours. A visual inspection was completed initially, at 30 minutes, one hour, and then every hour thereafter to 6 hours. The visual inspection included the following:

1. Verify mechanical integrity of ports
2. Verify mechanical integrity of the heat exchanger housing
3. Verify mechanical integrity of all seals and bonds

Any leakage was recorded.

### 17.2 Testing Results and Conclusion

The t = 0 (unaged) and t = 2 YR AA sample populations of the Affinity NT Oxygenator with Balance Biosurface built with Spiderman hollow microporous polypropylene fiber produced by Medtronic Danvers satisfied the water side integrity design verification performance requirements. Table 18-16 provides a summary of the results and conclusions.

Table 18-16: Water Side Integrity Testing

| Design Input Requirement (ID) | Acceptance Criteria                                               | %C / %R (Data Type) | Population  | n  | Result            | PASS/FAIL |
|-------------------------------|-------------------------------------------------------------------|---------------------|-------------|----|-------------------|-----------|
| Water Side Integrity (D1020)  | All devices must pass with no leaks over the six-hour time period | 95/90 (attribute)   | Unaged      | 29 | 29 of 29 no leaks | PASS      |
|                               |                                                                   |                     | t = 2 YR AA | 29 | 29 of 29 no leaks | PASS      |

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## 18 Blood and Gas Side Integrity and Port Break

### 18.1 Background

#### Pressure Integrity

The oxygenator is under constant fluid pressure and is responsible for heating and cooling the blood and transferring oxygen and carbon dioxide. There is pressure exerted on the blood-side of the device from the blood pump and patient. The blood-side must not leak for the oxygenator to function properly. As a result, a maximum rated operating pressure is specified for the bloodside. The integrity of the blood path is verified by subjecting the blood path to a hydraulic pressure more than the specified maximum

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recommended operating pressure for 6 hours, and then observing for leakage.

### Port Break

The device port strength must be greater than the force the ports may experience during foreseeable clinical use.

The intent of this test is to evaluate port break and blood side integrity performance after the gas path of the devices has been exposed to 1.3% Isoflurane and 2.6% Sevoflurane.

## 18.2 Testing Summary

### Blood Side Integrity

Deionized water was warmed to  $37 \pm 1^\circ\text{C}$  in a water bath. The pump was calibrated to a flow rate of  $3.0 \pm 0.2 \text{ L/min}$ . The blood path integrity circuit was constructed. Flow was initiated in the circuit at a flow rate of 3 LPM and continued at the required pressure for 6 hours. The circuit was pressurized using restrictor clamps to achieve the required pressure of  $23 \pm 1 \text{ psi}$  in the circuit. The pressure was confirmed at the inlet of the first device and the last device to ensure all are appropriately pressurized. A visual inspection was completed initially, at 30 minutes, one hour, and then every hour thereafter to 6 hours. The visual inspection included the following:

1. Verify mechanical integrity of ports and connectors.
2. Verify mechanical integrity of the housings.
3. Verify mechanical integrity of all bonds and seals.
4. Verify mechanical integrity of visible internal structures.

Any leakage was recorded. The test devices must not leak over the six-hour test period.

### Gas Side Integrity

The test circuit was constructed per Figure 18-6. A compressed air line was attached to the gas inlet port, and the gas outlet port was clamped. Air was slowly turned on to pressurize the device. The pressure was held for 10 minutes. If any leakage from the sealed gas outlet vents or any of the tubing stubs was observed, those leaks were repaired, and the test was repeated. A visual inspection was completed periodically, which includes the following:

1. Verify mechanical integrity of the ventilating gas cap
2. Verify mechanical integrity of the gas cap to oxygenator bond

Any leakage was recorded.  
Note: The same population of devices was used for gas and blood side pressure integrity testing. Gas side integrity testing was performed first.

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Table 18-17: Blood and Gas Side Integrity and Port Break Testing

| Requirement ID                   | Acceptance Criteria                     | %C / %R (data type) | Population   | n  | Result                                                    | PASS/FAIL |
|----------------------------------|-----------------------------------------|---------------------|--------------|----|-----------------------------------------------------------|-----------|
| Blood Path Integrity (DI019)     | No leaks over the six-hour test period  | 95/90 (attribute)   | T = 0        | 29 | No leaks                                                  | PASS      |
|                                  |                                         |                     | T = 2 YR AA  | 29 | 8 mandrel leaks <sup>1</sup>                              | FAIL      |
|                                  |                                         |                     | T = 18 MO AA | 29 | No Leaks                                                  | PASS      |
| Gas Path Integrity (DI021)       | No leaks over the 10-minute test period | 95/90 (attribute)   | T = 0        | 29 | No leaks                                                  | PASS      |
|                                  |                                         |                     | T = 2 YR AA  | 29 | No leaks                                                  | PASS      |
| Blood Port Break (DI026)         | Blood Port Break LTL ≥ 22.6 lbf         | 95/90 (variable)    | Unaged       | 11 | Blood Inlet: LTL = 217 lbf<br>Blood Outlet: LTL = 137 lbf | PASS      |
|                                  |                                         |                     | T = 2 YR AA  | 11 | Blood Inlet: LTL = 186 lbf<br>Blood Outlet: LTL = 142 lbf | PASS      |
| Gas Inlet Port Break (DI028)     | Gas Port Break LTL ≥ 22.6 lbf           | 95/90 (variable)    | Unaged       | 11 | LTL = 132 lbf                                             | PASS      |
|                                  |                                         |                     | T = 2 YR AA  | 11 | LTL = 108 lbf                                             | PASS      |
| Recirculation Port Break (DI029) | Recirc Port Break LTL ≥ 22.6 lbf        | 95/90 (variable)    | Unaged       | 11 | LTL = 122 lbf                                             | PASS      |
|                                  |                                         |                     | T = 2 YR AA  | 11 | LTL = 118 lbf                                             | PASS      |

<sup>1</sup> Eight mandrel leaks were observed in the two-year accelerated aged test population. After the mandrel joints of these devices were sealed per the work instructions (Attachment 12 of Attachment 18-20) and re-tested, no additional leaks were observed in all 29 samples over the six-hour duration of the test. See Section 18.18.1, Deviations and Test Failures for root cause investigation and re-test summary.

## 19 Plasma Breakthrough

### 19.1 Testing Summary

The device must be able to function throughout a case without blood plasma entering the hollow fiber. If this occurs gas cannot pass through the fiber and be transferred into (oxygen) and out of (carbon dioxide) the blood. Plasma resistance is achieved through the size of the pores, the hydrophobic nature of the fiber material, and any coating applied to the fiber. Individual fibers can fail if lipids penetrate and adhere to the hollow fiber pores allowing plasma to enter the inside of the hollow fiber. Therefore, plasma resistance is tested by flowing a lipid solution through the oxygenator and observing for plasma breakthrough in the fibers. Plasma breakthrough tests for blood plasma breaking through the hollow fibers of the oxygenator's hollow fiber bundle for the extent of the six-hour indication of the device. This test is intended to assess plasma breakthrough performance after the gas path of the

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devices have been exposed to 1.3% Isoflurane and 2.6% Sevoflurane.

Lipid solution was circulated through the device at a maximum rated flow for six hours. The device was observed during the test to determine if breakthrough had occurred. Greater than 60 cumulative plasma breakthrough fiber leaks at the end of the six-hour test was considered plasma breakthrough. The test circuit was set up.

## 19.2 Testing Results and Conclusion

The t = 0 (unaged) and t = 2 YR AA sample populations of the Affinity NT Oxygenator with Balance Biosurface built with Spiderman hollow microporous polypropylene fiber produced by Medtronic Danvers satisfied the plasma breakthrough design verification performance requirements after 6 hours of gas path exposure to 1.3% Isoflurane and 2.6% Sevoflurane. Table 18-19 provides a summary of the results and conclusions.

Table 18-19: Plasma Breakthrough Testing

| Design Input Requirement (ID) | Specification Acceptance Criteria                                                                                  | %Confidence / %Reliability (Data Type) | Population  | Device                  | n              | Result (μ ± σ)           | PASS/FAIL |
|-------------------------------|--------------------------------------------------------------------------------------------------------------------|----------------------------------------|-------------|-------------------------|----------------|--------------------------|-----------|
| DI018                         | Devices shall not exhibit plasma breakthrough over a six-hour test period under blood and gas flow rates of 7 LPM. | 90/90 (attribute)                      | t = 0       | Balance NT (Test)       | 22             | 22 of 22 no breakthrough | PASS      |
|                               |                                                                                                                    |                                        |             | Trillium NT (Predicate) | 4              | N/A = predicate device   |           |
|                               |                                                                                                                    |                                        | t = 2 YR AA | Balance NT (Test)       | 22             | 22 of 22 no breakthrough | PASS      |
|                               |                                                                                                                    |                                        |             | Trillium NT (Predicate) | 3 <sup>1</sup> | N/A = predicate device   |           |

<sup>1</sup> On the third test day, testing of the control device C3 (serial number 5116124480) had to be stopped 1 hour and 15 minutes into the test due to a mandrel leak. See Section 18.19.3 Deviations below for details.

## 20 Luer Torque

### 20.1 Testing Summary

The intent of this test to measure the device luer thread strength, which must be greater than the force the threads may experience during foreseeable clinical use.

A rack and pinion fixture were secured to Instron platform, such that the straight-toothed rack is level. The test device is then placed in a cradle, and a polycarbonate male luer fitting is attached to the test port. The test port is aligned and inserted into the rack and pinion fixture on the Instron.

To apply a torque to the luer port under test, the Instron moved the straight-toothed rack at a constant displacement, which was converted into rotational motion of the polycarbonate male luer

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fitting attached to the test port. The maximum torque applied to the test port was recorded.

## 20.2 Testing Results and Conclusion

The  $t = 0$  (unaged) and  $t = 2$  YR AA sample populations of the Affinity NT Oxygenator with Balance Biosurface built with Spiderman hollow microporous polypropylene fiber produced by Medtronic Danvers satisfied the luer torque design verification performance requirements. Table 18-20 provides a summary of the results and conclusions.

Table 18-20: Luer Torque Testing

| Design Input Requirement (ID) | Specification Acceptance Criteria               | %Confidence / %Reliability (Data Type) | Population    | Device                     | n  | Result ( $\mu \pm \sigma$ ) | PASS/FAIL              |
|-------------------------------|-------------------------------------------------|----------------------------------------|---------------|----------------------------|----|-----------------------------|------------------------|
| Luer Torque (D1030)           | The luer port torque shall be $\geq 5.5$ in-lbf | 95/90 (variable)                       | Unaged        | Balance NT - Access Port   | 11 | $20.3 \pm 1.7$ in-lbf       | LTL = 18.9 in-lbf PASS |
|                               |                                                 |                                        |               | Balance NT - Sampling Port | 11 | $27.4 \pm 3.0$ in-lbf       | LTL = 23.2 in-lbf PASS |
|                               |                                                 |                                        | $t = 2$ YR AA | Balance NT - Access Port   | 11 | $22.9 \pm 3.1$ in-lbf       | LTL = 15.7 in-lbf PASS |
|                               |                                                 |                                        |               | Balance NT - Sampling Port | 11 | $24.6 \pm 4.4$ in-lbf       | LTL = 18.0 in-lbf PASS |

## 21 Conclusion

The test results documented in sections above and summarized within Table 18-21 for the subject devices, Balance coated Affinity NT Oxygenator (standalone – 511B and combo w/uncoated CVR – 541B), demonstrate that the performance is comparable to the predicate Affinity NT Oxygenator devices (511T and 541T) currently in market.

Table 18-21: Testing Results

| Test Description                              | T=0 Result | T=2yrAA Result    |
|-----------------------------------------------|------------|-------------------|
| Balance Coverage                              | PASS       | PASS              |
| Balance Leaching                              | PASS       | PASS              |
| Initial Gas Transfer & Sampling Port Accuracy | PASS       | PASS              |
| Time Dependent Gas Transfer                   | PASS       | PASS              |
| Prime Volume                                  | PASS       | PASS              |
| Blood Trauma (maximum flow)                   | PASS       | PASS              |
| Blood Trauma (minimum flow)                   | PASS       | PASS              |
| Plasma Breakthrough                           | PASS       | PASS              |
| Pressure Integrity (blood, gas)               | PASS       | PASS <sup>1</sup> |
| Pressure Integrity (water)                    | PASS       | PASS              |
| Burst (blood, gas)                            | PASS       | PASS              |

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| Description                      | T=0 Result | T=2yrAA Result |
|----------------------------------|------------|----------------|
| Balance Coverage                 | PASS       | PASS           |
| Burst (water)                    | PASS       | PASS           |
| Port Break                       | PASS       | PASS           |
| Tubing Pull                      | PASS       | PASS           |
| TMA Accuracy                     | PASS       | PASS           |
| Dust Cap Removal                 | PASS       | PASS           |
| Particulate                      | PASS       | PASS           |
| * Tested at 18-months shelf life |            |                |

## Overall Conclusion

In conclusion, the information included in this submission demonstrates that the subject devices, Affinity NT Oxygenator with Balance Biosurface (standalone - 511B and combo configuration integrated with uncoated CVR - 541B) built with Medtronic manufactured fiber, are substantially equivalent to the legally marketed predicate devices (511T and 541T).

## Biocompatibility

| Description                          | Standard     | Result |
|--------------------------------------|--------------|--------|
| Cytotoxicity                         | ISO 10993-5  | PASS   |
| Sensitization                        | ISO 10993-10 | PASS   |
| Irritation / Intracutaneous Toxicity | ISO 10993-10 | PASS   |
| Systemic toxicity                    | ISO 10993-11 | PASS   |
| Hemocompatibility                    | ISO 10993-4  | PASS   |

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
*Pressure Decay Test (31058)*

*Leak test*

*Place Fiber bundle sub-assembly on leak tester nest*

*Place FBA in to Leak tester Nest*

*Fiber Leak Test Units*

*Inspect both sides of every fiber bundle during test, looking for water droplets to form*

*Fiber Leak Test (31187)*

*Drying Fiber Bundle after Leak Test (31204)*

*Inspect for moisture, water spots, damaged ports; retry if moisture present. Final Outer term should at least 30 deg C*

*Dry Units*

*Place FBA in to Dryer Nest*

**LEGEND**

*Manufacturing Processes performed at Brooklyn Park, MN facility*

*Manufacturing Processes performed at Brooklyn Park, MN facility & Tijuana, Mexico*

*Sterilization at Vention Sterilization Service, Grand Rapids, MI & Midwest Sterilization, Jackson, Mo*

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
#### HEAT EXCHANGE SUB-ASSEMBLY

Place lower housing on nest, then position below in it. Then place upper housing on the subsequent nest

Dispense UV adhesive on lower housing groove and upper housing installed

Cure

HEX side seam Bond (31171)

Water Diverter potting (31168)

Place Hex Assembly in oven for curing

Place Hex assembly on nest oriented for water diverter potting

Close End potting (31155)

Place Hex Assembly in oven for curing

Place Hex assembly on nest oriented for Closed End potting

Hex Leak Test (31171)

Leak test

Place Hex Assembly on leak tester nest

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

**Main Bond (31174)**  
 Cure  
 Dispense UV Adhesive  
 Load HFO assembly in station  
 Trillium coating  
 2x Rinse cycle  
 Polymer B coat cycle  
 2x Rinse cycle  
 Polymer A coat cycle  
 Prepare Polymer A & B solution  
 Place counted units in dryer  
 Dry units  
 Drying  
 Sterile Packaging  
 Seal bags  
 Integrate with CVR  
 Sterilization

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
### AFFINITY NT OXYGENATOR (HFO) MANUFACTURING FLOWCHART

# **FIBER BUNDLE SUB-ASSEMBLY**

## **FIBER BUNDLE SUB-ASSEMBLY**

*Mandrel Preparation*

*Install Bobin*

*Set winding*

*Wind*

*Assemble Fiber bundle (31184)*

*Fiber bundle Potting (31185)*

*De-cap Fiber bundle*

*Remove units and Put small end*

*Place potting cap feature in centrifuge (Large End), add Urethane leveler and cap*

*Set Centrifuge*

*Assemble Potting fixture: small End Cap, liners, Oxy Out case, bundle and large end Cap*

*Prep. Potting Cap features*

*Pre-heat Fiber Bundle*

*Remove Small End Liner*

*Guideline Excess material*

*Guillotined Fiber Bundle (31186)*

*Gas Cap Bond (31202)*

*Cure*

*Install Gas Cap*

*Dispense UV adhesive in outer case gas cap groove*

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

Manufacturing performed at Brooklyn Park, MN

Manufacturing process performed at Brooklyn Park, MN & Tijuana, Mexico

Pressure Decay Test (31058)  
Leak test  
Place Fiber bundle sub-assembly on leak tester nest

Place FBA in to Leak tester Nest

Fiber Leak Test Units

Inspect both sides of every fiber bundle during test, looking for water droplets to form

Fiber Leak Test (31187)

Drying Fiber Bundle after Leak Test (31204)

Inspect for moisture, water spots, damaged ports; retry if moisture present. Final Outer temp should at least 30 deg C

Dry Units

Place FBA in to Dryer Nest

Manufacturing performed at Brooklyn Park, MN

Manufacturing process performed at Brooklyn Park, MN & Tijuana, Mexico

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
**Pressure Decay Test (31058)**

**Leak test**

**Place Fiber bundle sub-assembly on leak tester nest**

**Place FBA in to Leak tester Nest**

**Fiber Leak Test Units**

**Inspect both sides of every fiber bundle during test, looking for water droplets to form**

**Fiber Leak Test (31187)**

**Drying Fiber Bundle after Leak Test (31204)**

**Inspect for moisture, water spots, damaged ports; retry if moisture present. Final Outer term should at least 30 deg C**

**Dry Units**

**Place FBA in to Dryer Nest**

**LEGEND**

**Manufacturing Processes performed at Brooklyn Park, MN facility**

**Manufacturing Processes performed at Brooklyn Park, MN facility & Tijuana, Mexico**

**Sterilization at Vention Sterilization Service, Grand Rapids, MI & Midwest Sterilization, Jackson, Mo**

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
Main Bond (31174)  
 Cure  
 Dispense UV Adhesive  
 Load HFO assembly in station  
 Balance coating  
 2x Rinse cycle  
 Polymer B coat cycle  
 2x Rinse cycle  
 Polymer A coat cycle  
 Prepare Polymer A & B situation  
 Place counted units in dryer  
 Dry units  
 Drying  
 Sterile Packaging  
 Seal Bags  
 Integrate with CVR  
 Sterilization

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Affinity NT oxygenator with Trillium biosurface contains heparin. The Trillium Biopassive surface is a hydrophilic polymer coating with heparin that increase hemocompatibility of the device. Background: Trillium Biosurface are biocompatible coatings applied to devices marketed by Medtronic Cardiovascular RST. These coatings contain heparin derived from porcine intestinal mucosa, which is obtained from manufacturers of heparin that adhere to the USP heparin sodium monograph standard. Risk Management of heparin in Trillium coating: As described in ISO 22442-1:2007, hazards could arise from  
 a) contamination by bacteria, moulds or yeasts;

**Pharma product  
contained in medical  
device (if applicable)**

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b) contamination by viruses;

c) contamination by agents causing Transmissible Spongiform Encephalopathies (TSE);

d) material responsible for undesired pyrogenic, immunologic or toxicological reactions.

Hazards arising from b) and c) could potentially lead to patient infections that could be irreversible and, therefore, be considered high severity. To mitigate this risk, the USP heparin monograph standard requires manufacturers of heparin to certify that appropriate measures have been taken to ensure that the heparin being supplied is free of contamination by viruses and TSE agents. The certificate of analysis provided by the heparin sodium manufacturer includes information on levels of bacterial endotoxins. Additional risk mitigation is achieved through sterilization of the coated medical device and testing to comply with ISO 10993. Based on these, the residual risk for the use of heparin sodium in the Trillium coating is deemed acceptable.


**Heparin Sodium**

Chemical Abstract Service Number: 9041-08-1

Structure (Heparin Sodium, Smithfield BioScience Inc.)

Heparin Sodium is a substance comprising the sodium salt of sulfated glycosaminoglycans isolated from mammalian tissue. On complete hydrolysis it liberates D-glucosamine, D-glucuronic acid, L-iduronic acid, acetic acid and sulfuric acid. It is polydisperse in chain length and heterogeneous in degree and type of saturation. Much of the structure contains repeating trisulfated disaccharides, disulfated disaccharides, and disaccharides lacking a sulfate on C-6 of the glucosamine residue.

The general chemical structure of Heparin is shown below:

**General Structure of Heparin**

The anticoagulant properties of heparin are primarily mediated by a specific pentasaccharide structure with high affinity for the serine protease inhibitor antithrombin (AT), a structure of which comprises no more than about 30% of the total heparin oligosaccharide composition. The high-affinity AT-binding pentasaccharide is shown below:

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
#### High-Affinity AT-binding Polysaccharide

**Мolecular Formula:** The molecular formula is  $(C_{12}H_{16}NS_2Na_3)_n$ . Heparin sodium does not contain unique stereochemistry specific to its activity.

**Мolecular Weight:** Heparin Sodium is the sodium salt of sulfated glycosaminoglycans present as a mixture of heterogeneous molecules varying in molecular weight of polyanion chains, and therefore has range of molecular weights ( $M_w$ ). Elemental analysis of typical batches of heparin Sodium results in 24-25% carbon, 3.4-3.6% hydrogen, 1.9-2.3% nitrogen, and 9.5-11.0% sulfur. Activity measurements made with an ion-specific electrode show that only 66% of the sodium ions in heparin sodium are bound.

The molecular weight range of Heparin Sodium USP 37 defined as follows:

- $M_w$  above 24 kDa is not more than (NMT) 20%
- Average  $M_w$  is between 15 kDa and 19 kDa.
- The ratio of the percentage  $M_w$  8-16 to  $M_w$  16-24 is Not Less Than (NLT) 1.0.

**General Properties (Heparin Sodium, Smithfield BioScience Inc.)**

**Physicochemical Properties:**

Medtronic uses heparin sodium which meets the compendial requirements of the United States Pharmacopeia (USP) and the European Pharmacopocia (Ph. Eur.). Heparin sodium is a sodium salt powder with a white to off-white appearance. Heparin sodium is soluble in water up to 60% by mass, but essentially insoluble in methanol, ethanol or acetone. Heparin sodium is a clear liquid in a 10% solution.

The pH of a 1% solution of Heparin Sodium is 5.5 – 7.5.


For Heparin Sodium, World Health organization (WHO) International Standard, 6th International Standard for Unfractionated heparin, NIBSC code: 07/328, 2145 IU/ampoule, is being employed by Smithfield BioScience Inc. in the determination of potency in IU per mg.

The physiochemical properties of Heparin Sodium are provided in report, held by the manufacturer, Smithfield BioScience Inc. and filed with MHRA.

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
Manufacturer(s) (Heparin Sodium, Smithfield BioScience Inc.)  
 The manufacturing information is as follows:  
 Drug substance, heparin sodium:  
 Smithfield BioScience Inc.  
 12150 Best Place  
 Cincinnati, Ohio 45241  
 United States of America

Description of Manufacturing Process and Process Controls (Heparin Sodium, Smithfield BioScience Inc.)  
 The heparin sodium manufactured by Smithfield BioScience Inc. is sourced from both the USA and China as indicated by the figure below.

Manufacturing Flowchart Including Heparin Sources

Description and Composition of Drug Product (Trillium Biosurface)

Trillium Biosurface (Trillium) is a heparin-based surface-modification technology that mimics the heparan sulfate naturally found on the vascular endothelium lining of the circulatory system. The Trillium coating is applied to the artificial surface of various cardiopulmonary medical devices. A complete list of medical devices and respective coated artificial surfaces.

Trillium is comprised of the following components:

- Polymer A

The priming layer, Polymer A, consists of polyethyleneimine (PEI). The solid content of polymer A received by Medtronic is 5.8% - 7.8% as documented on the Certificate of Analysis. Polymer A is diluted with RODI water prior to coating of the devices to a specified dilution. The priming layer is deposited on the device surface by perfusing the diluted Polymer A through the device.

- Polymer B (with Heparin)

The functional layer, Polymer B, consists of Polyethylene Oxide (PEO) chains on to which heparin is covalently attached. The PEO chain has a high degree of sulfonation (negative charge). The solid content of the Polymer B received by Medtronic is 1%- 2% as

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
(АТIII). Heparin binds to АТIII and initiates a change in АТIII that leads to the inactivation of Thrombin and FXa (an additional component of coagulation). This effect has been demonstrated as well in terminally sulfated PEO chains to the degree of 14% of the anticoagulant capacity of heparin in in-vitro studies.

#### Trillium™ Biosurface

Depiction of Trillium Coating structure on Medical Devices

#### Functional Layer


The functional layer of Trillium coating is referred to as Polymer B. At BioInteractions, Ltd., the heparin sodium drug substance purchased from Smithfield BioScience Inc. is reacted with activated Polyethylene Glycol (6) Methacrylate (PEGM) and pH adjusted to produce the Polymer B.

#### Priming Layer

The priming layer for Trillium coating referred to as Polymer A consists of a linear polyethyleneimine (PEI) modified with additional alkyl chemistry to render the polymer's hydrophobic character and promote adhesion to similar hydrophobic substrates e.g., disposable plastic devices; the amide component confers some cationic character to the alkyl modified PEI that allows it to bond to other charged surfaces such as metals and serves as the anchor for the functional top layer (Polymer B).

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**Affinity NT oxygenator with Balance biosurface, including Warnings**

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■ The fluid pathway is sterile and nonpyrogenic. Inspect each package and device before use. Do not use the device if the package is opened or damaged, the device is damaged, or the protective caps are not in place.

■ Continually monitor the extracorporeal circuit for air bubbles or leaks. Do not use the device if air bubbles or leaks are observed during priming or operation as they may result in air embolism to the patient or fluid loss.

■ Do not operate the oxygenator and CVR outside the recommended parameter ranges.

■ Once an oxygenator is primed with blood, adequate anticoagulation should be maintained per institution CPB protocol and the blood pathway should be constantly recirculated at a low flow rate.

■ After initiation of bypass, blood flow should circulate through the oxygenator and CVR at all times within the recommended blood flow range during CPB, except during emergency oxygenator replacement.

■ Ensure that unused ports remain capped to prevent contamination and that the caps are tightened to avoid leakage.

■ Clear all gas emboli from the extracorporeal circuit before initiating bypass. Gas emboli are hazardous to the patient.

■ In any extracorporeal circuit, the use of a prebypass and arterial line filter is recommended.

■ Monitoring of circuit pressure is recommended.

■ The blood phase pressure must be greater than the gas phase pressure at all times in the oxygenator.

■ Do not exceed blood pressure of 100 kPa (750 mmHg) in the blood pathway of the oxygenator.

■ Do not exceed heat exchanger water phase pressure of 206.8 kPa (30 psi).

■ Heat exchanger water temperature shall not exceed 42°C (107°F).

■ Follow institution CPB protocol when using disinfectants in the heater/cooler during bypass. The integrity of the water path has been verified with hydrogen peroxide (330 ppm). Contact

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- May be used for information regarding the use of additional disinfectants.
- Band all tubing connections for added protection against high fluid pressures.
  - A replacement oxygenator should be readily available during CPB.
  - Do not obstruct gas outlet vents to prevent the gas side pressures from exceeding blood side pressures.
  - Gas transfer rates may change over time, and adjustments of FiO2 and gas flow rates may be needed to achieve desired gas transfer performance.
  - The use of a blood level sensing mechanism during operation of this device is recommended.
  - Ensure that the CVR outlet is always positioned above the highest point in the membrane compartment of the oxygenator.
  - Do not use a vacuum source in excess of 150 mmHg suction through this device.
  - Do not occlude or obstruct the vent port of the CVR during operation (except when using VAVD).
  - Only persons thoroughly trained in CPB procedures should use this device. Operation of each device requires constant supervision by qualified personnel for patient safety. CPB must be carefully and constantly monitored.
  - This device was designed for single patient use only. Do not reuse, reprocess, or resterilize any component of the system. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device, which could result in patient injury, illness, or death.
  - Each device is sterilized using ethylene oxide.
  - Attach tubing in such a manner as to prevent kinks or restrictions that may alter blood, water, or gas flow.
  - Do not allow alcohol, alcohol-based fluids, anesthetic liquids (such as isoflurane), or corrosive solvents (such as acetone) to come into contact with the device as they may jeopardize the structural integrity.

#### Precautions

- Use aseptic technique in all procedures.
- Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures. Maintain adequate anticoagulation per institution CPB protocol.

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
1. Cardiotomy/venous reservoir
2. Oxygenator
3. Arterial filter
4. Prebypass filter

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**Warning:** Ensure that aseptic technique is used during all stages of setup and use of this system.

**Warning:** Before removing the devices from the packaging, inspect the packaging and product for damage.

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
**Note:** The Affinity NT oxygenator and CVR may be separated and utilized with separate holders. To separate the oxygenator from the CVR, turn the oxygenator approximately 1/8 of a turn in a counterclockwise direction and place in their respective holders.

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
5. Connect all blood and gas lines to their appropriate connections. Do not occlude tubing in pump raceway.
6. Clamp the tubing at the oxygenator inlet, oxygenator arterial outlet, and oxygenator recirculation line outlet.


**Note:** A CO<sub>2</sub> flush of the oxygenator is not required, but is recommended before priming.

**Priming (Refer to figure above)**

1. Fill the CVR with 1.5 to 2 L of priming solution.
2. Remove all bubbles in the section of tubing from the CVR outlet to the clamp by the oxygenator inlet.
3. Remove the clamp from the oxygenator inlet and allow solution to gravity prime the heat exchanger and fiber bundle of the oxygenator.

**Note:** Albumin has been reported to reduce the rare instance of transient high inlet pressure phenomenon in membrane oxygenators. Therefore, albumin in the prime solution should be considered.

**Recirculation (Refer to Figure above)**

1. After placing tubing into a properly occluded pump raceway, remove clamp from recirculation outlet per Figure above to allow recirculation between the oxygenator and venous reservoir.

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
2. Begin recirculation flow at 500 ml/minute and gradually increase flow to 5 liters/minute.

3. Ensure that oxygenator is bubble free.

4. Gradually reduce flow, stop pump and clamp recirculation line.

Priming A/V Circuit (Refer to Figure above)

1. Add additional prime solution to fill the rest of the perfusion circuit if needed. Unclamp the arterial outlet; slowly pump prime solution into the remaining circuit. Purge air from the arterial sample port and access port at this time.
2. The arterial filter should now be primed per manufacturer's instructions.
3. Inspect the entire circuit and all components for leaks.
4. Upon completion of priming and debubbling, gradually reduce blood flow and stop the pump, close all purge lines, and clamp the arterial, venous, and recirculation lines.
5. Ensure that the entire extracorporeal circuit is free of bubbles before initiating bypass.

Initiation of Bypass

Warning: The blood phase pressure must be greater than gas phase pressure at all times.

Note: Priming solution may be pre-warmed through heat exchanger prior to initiation of bypass.

1. Check for adequate levels of anticoagulation prior to initiation of bypass.
2. Remove arterial and venous clamps and gradually increase blood flow. Then initiate gas to blood flow ratio of 1:1 with a minimum FiO<sub>2</sub> of 80%.
3. Initiate water flow through heat exchanger.

Warning: Do not allow air to enter the circuit.

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6. To remove the oxygenator from the CVR, turn the oxygenator approximately 1/8 of a turn in a counterclockwise direction.
7. Attach the replacement oxygenator to the CVR by inserting the positioning pegs located on the gas cap into the coupling slots on the CVR, and turn clockwise until the oxygenator locks into position.
8. Replace the heat exchanger support base onto the holder by sliding the base back up the mast until the heat exchanger rests securely in the horizontal channel of the base. Secure the support base by inserting the release pin of the support base into the holder mast.
9. Reconnect the recirculation, oxygenator inlet, and outlet tubing. Band all connections.
10. Connect water lines, unclamp, turn on water source, and inspect for water leaks.
11. Connect oxygen line.
12. Remove clamps from recirculation and oxygenator inlet tubing. Warning: Do not unclamp arterial line or venous return line at this time.
13. With sufficient volume in the CVR, turn on blood pump and slowly fill the oxygenator.
14. Increase blood flow to 5 L/min through recirculation line.
15. Stop pump and clamp recirculation line.
16. Ensure that the entire system is free of leaks and gas bubbles.
17. Remove all clamps from venous and arterial lines and reinitiate bypass.

**Affinity NT oxygenator with Balance biosurface, including System setup**

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Warning: Before removing the devices from the packaging, inspect the packaging and product for damage.

If the package or product is damaged, do not use as sterility of the device may have been compromised and/or performance may be affected.

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Note: The Affinity NT oxygenator and CVR may be separated and utilized with separate holders. To separate the oxygenator from the CVR, turn the oxygenator approximately 1/8 of a turn in a counterclockwise direction and place in their respective holders.

3. Remove the obturator cap from the positive/negative pressure relief valve.

4. Connect the 1.3 cm (1/2 in) water lines to the inlet and outlet ports on the heat exchanger. Start the heat exchanger water flow, and check for leaks from the water compartment to the blood compartment. Discard

the oxygenator if water is present in blood compartment, and set up the standby oxygenator.

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6. Clamp the tubing at the oxygenator inlet, oxygenator arterial outlet, and oxygenator recirculation line outlet.

Note: A CO<sub>2</sub> flush of the oxygenator is not required, but is recommended before priming.

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9. Connect the quick-connect 1.3 cm (1/2 in) water lines to the inlet and outlet ports on the oxygenator. Start the water flow and check for leaks from the water compartment to the blood compartment before priming.

Caution: Do not use the oxygenator if water is present in the blood compartment.

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replace the heat exchanger support base onto the holder by sliding the base back up the mast until the heat exchanger rests securely in the horizontal channel of the base. Secure the support base by inserting the release pin of the support base into the holder mast.

9. Reconnect the recirculation, oxygenator inlet, and outlet tubing. Band all connections.
  10. Connect the water lines, unclamp, turn on the water source, and inspect for water leaks.
  11. Connect the oxygen line.
  12. Remove the clamps from the recirculation and oxygenator inlet tubing.
- Warning:** Do not unclamp the arterial line or the venous line at this time.
13. With sufficient volume in the CVR, turn on the blood pump and slowly fill the oxygenator.
  14. Increase blood flow through the recirculation line.
  15. Stop the pump and clamp the recirculation line.
  16. Ensure that the entire system is free of leaks and gas bubbles.

#### VAVD instructions for use

1. Follow previously listed instructions for use for setup with the following modifications as listed here in steps 2 through 3.
  2. Attach a calibrated means of monitoring negative pressures to either the venous line or the CVR.
  3. Follow previously listed instructions for use for priming with the following additions as listed here in steps 5 through 9.
  4. Replace the vented white luer caps on the manifold sampling line ports with nonvented luer caps. Place plugs ("dead-enders") on all unused 0.6 cm (1/4 in) and 1.0 cm (3/8 in) ports.
  5. Prime, debubble, and complete preparation of the circuit as per standard procedure. Prime and zero the venous line pressure monitor.
  6. Before instituting bypass, prepare the vacuum portion of the setup. Attach a vented vacuum line to the 0.6 cm (1/4 in) CVR vent port and to the vacuum regulator with vapor trap.
- Warning:** Clamp the quick prime line before application of vacuum.
- Warning:** Do not allow the vapor trap to become completely filled during use. This may allow fluid to enter the vacuum controller, or may prevent the CVR from being vented to atmosphere when no vacuum is applied.

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|                                  |                                                                                                                                                           |
|----------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| Directive 2001/83/EC             | Directive 2001/83/EC of the European Parliament and of the Council of 6 November 2001 on the Community code relating to medicinal products for human use. |
| European Pharmacopeia (Ph. Eur.) | European Pharmacopeia 7th Edition (v7.8) sections 3A, 3B, 3C                                                                                              |

**Affinity NT oxygenator with Balance biosurface, including**

| Standard Number | Description                                                                                                  |
|-----------------|--------------------------------------------------------------------------------------------------------------|
| EN ISO 13485    | Medical Devices - Quality Management Systems - Requirements for regulatory purposes                          |
| ISO 15675       | Cardiovascular implants and artificial organs - Cardiopulmonary bypass systems - Arterial blood line filters |
| EN ISO 14971    | Application of Risk Management to Medical Devices                                                            |
| ISO 7199        | Cardiovascular implants and artificial organs - Blood gas exchanges (oxygenators)                            |
| EN 1041         | Information Supplied by the Manufacturer with Medical Devices                                                |
| ISO 15223-1     | Medical Devices - Symbols to be Used with Medical Device Labels, Labeling and Information to be Supplied     |
| EN ISO 10993-1  | Biological Evaluation of Medical Devices - Part 1: Evaluation and Testing within a risk management process   |
| EN ISO 10993-2  | Biological evaluation of medical devices - Part 2: Animal welfare requirements                               |

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© 2018 Medtronic