ZEISS

Carl Zeiss Meditec AG 73446 Oberkochen

Roszdrazvnadzor

Carl Zeiss Meditec AG

Oberkochen Site
Rudolf-Eber-Str. 11
73447 Oberkochen

Phone: +49 (0) 7364/ 20 8613

e-mail: anke.seitz@zeiss.com

Division/Dept.: MED-RAO
Your contact: Anke Seitz

Your ref.:

Yours of:

Our ref.:

Date: 04.05.2021

Instruction for Use (IFU)

Dear Madam or Sir,

Please find enclosed the English and Russian Instruction for Use (IFU) for KINEVO 900.
G_30_2094_RU_V1_2_KINEVO_900 + addendum

G_30_2094_v1.1_en_KINEVO_900_SW1.8 + addendum

Best regards

Carl Zeiss Meditec AG
i.V. Anke Seitz
Regulatory Affairs Manager

Regulatory
Affairs

Dr. Jürgen Berger
Regulatory Affairs Manager

Manufacturing Plant Oberkochen
Rudolf-Eber-Strasse 11
73447 Oberkochen, Germany

Address of Record:
Goeschwtitzer Strasse 51 - 52
07745 Jena, Germany

Address for Delivery:
Carl Zeiss Meditec AG
Rudolf-Eber-Strasse 11
73447 Oberkochen, Germany

Banks:
Deutsche Bank Jena
Account: 624536900 (BLZ 820 700 00)
IBAN: DE90 8207 0000 0624 5369 00
BIC/ SWIFT: DEUT DE 8EXXX

Commerzbank Jena
Account: 258072800 (BLZ 820 400 00)
IBAN: DE31 8204 0000 0258 0728 00
BIC/ SWIFT: COBADEFFXXX

Commercial Register:
Local Court Jena HRB 205623

VAT-ID No.: DE 811 922 737
WEEE-Reg.-No.: DE55298748

Chairman of the Supervisory Board:
Dr. Michael Kaschke

Board of Management:
Dr. Ludwin Monz (CEO)
Justus Felix Wehmer
Jan Willem de Cler

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

© 2020, Carl Zeiss Meditec AG, Jena

Brands

ZEISS KINEVO 900 is a trademark or registered trademark of Carl Zeiss Meditec AG or other companies of the ZEISS Group in Germany and/or other countries.

All other brands mentioned in this document are the property of their respective owners.

Table of Contents

1	Notes on the instructions for use	11
1.1	Product name	11
1.2	Scope of application	11
1.3	Purpose and storage of the documentation	11
1.4	Questions and comments	11
1.5	Conventions in this document	11
1.5.1	Conventions in all text areas	11
1.5.2	Conventions in a course of action	12
1.6	Other applicable documents	13
2	Safety notes	15
2.1	Target group	15
2.2	Area of use	15
2.2.1	Intended use	15
2.2.2	Indications	16
2.2.3	Normal use	16
2.2.4	Possible damage/injury from burns due to the high light intensity	18
2.2.5	Electromagnetic compatibility	21
2.3	Responsibilities and duties of the operator	21
2.3.1	Requirements for the network connection	23
2.3.2	Messages to manufacturer and authorities	25
2.4	Measures and duties of the operator	25
2.5	Liability and warranty	26
3	Device description	27
3.1	General	27
3.2	Innovative functionalities	27
3.2.1	Kinematic and robot-assisted positioning functions	27
3.2.2	QEVO handheld Micro-Inspection Tool (option)	28
3.2.3	Integrated 3D video system for observation without eyepieces (option)	28
3.2.4	Data management	29
3.2.5	Navigation functionality	29
3.2.6	Video recording functionality (option)	30
3.3	Device marking	30
3.3.1	Labeling on the microscope	31
3.3.2	Approval label / UDI label	32
3.3.3	Labeling on stand, Part 1	33
3.3.4	Labeling on stand, Part 2	34
3.3.5	Labeling on stand, Part 3	36
3.3.6	Labeling on stand base	38
3.3.7	Labeling on light source	39
3.3.8	Labeling on autofocus option	40
3.3.9	Labeling on packaging	41

3.3.10	Displaying radio frequency approvals on the monitor	42
3.3.11	Display Windows EULA on the monitor	43
3.4	Structure of the device	44
3.4.1	Device overview	44
3.4.2	Overview of connector panel	45
3.4.3	Overview of video connector panel	47
3.4.4	Overview of video connector panel (4K 2D/3D option)*	49
3.4.5	Overview of microscope	50
3.5	Control elements and displays	56
3.5.1	Screen keyboard	56
3.5.2	Display: Editing	56
3.5.3	Online help	57
3.5.4	Error messages on the monitor	58
3.5.5	Controls on the microscope	59
3.5.6	Controls of binocular tubes	60
3.5.7	Controls for widefield eyepieces	61
3.5.8	Controls of handgrips	62
3.5.9	14-function foot control panel	63
3.6	Software description	64
3.6.1	General	64
3.6.2	Interactive buttons	64
3.6.3	User interface of the KINEVO 900	66
3.6.4	QEVO user interface (option)	69
3.6.5	Status Information	70
3.6.6	"Patient data" menu	71
3.6.7	Importing patient data from the DICOM Worklist	73
3.6.8	"Edit video recording" menu	75
3.6.9	SmartLinks	77
3.6.10	Light	78
3.6.11	Tube	79
3.6.12	Focus	80
3.6.13	Zoom	81
3.6.14	Hand Grips	82
3.6.15	Foot Control Panel (FCP)	84
3.6.16	Rock Switch	85
3.6.17	Suspension system	86
3.6.18	Displays	87
3.6.19	Audio	88
3.6.20	Photo	89
3.6.21	Video Camera	90
3.6.22	Connections	93
3.6.23	Video recording	102
3.6.24	Fluorescence	103
3.6.25	MultiVision	106
3.6.26	Info	108
3.6.27	Extras	109
3.6.28	Service PC	110
4	Installation	115
4.1	Requirements	115
4.2	Attaching the tubes and eyepieces	116

4.3	Using an integrated 3D video system for observation without eyepieces (option).....	118
4.4	Attaching the documentation/co-observer facility	120
4.5	Mounting the mouth switch	121
4.6	Moving the device	122
4.7	Connecting the device to the power supply.....	124
4.8	Switching the device on	126
4.9	Initial device login	127
4.10	Connecting the wired foot control panel or rocker foot switch	128
4.11	Pairing the wireless foot control panel with the device	129
4.12	Connecting video systems and external monitors.....	132
4.12.1	Connecting video devices	134
4.12.2	Connection of the 4K 3D monitor	134
4.12.3	Connection of the 4K 2D monitor	138
4.12.4	Connection of the HD 3D monitor	142
4.12.5	Connection of the HD 2D monitor	145
4.12.6	Connecting the monitor of a third-party manufacturer	148
4.13	Connecting USB storage media	148
4.14	Configure network	149
4.14.1	Activating IT admin system rights	149
4.14.2	Activate network connection via LAN	150
4.14.3	Activate network connection via wireless LAN (option).....	151
4.14.4	Creating and activating the device as a hotspot.....	151
4.14.5	Activating the web interface in the device	152
4.14.6	Network drive	152
4.14.7	Configure the network connection for DICOM (option).	153
4.14.8	Changing the computer name of the device	153
4.15	Configuring the Service PC	154
4.15.1	Service	154
4.15.2	Settings	154
4.15.3	Log Files	156
4.16	Connecting the navigation system to the device	156
4.17	Mounting the laser micromanipulators.....	157
5	Daily startup.....	159
5.1	Safety during preparation.....	159
5.2	Setting the position of the handgrips	160
5.3	Positioning the device in the OR.....	161
5.4	Check that the system can move freely!	161
5.5	Configuring the device for applications	162
5.6	Drape position	162
5.7	Attaching a ZEISS SMARTDRAPE.....	163
5.8	Autobalance	165
5.9	Setting the microscope.....	166

5.10	Connect QEVO Mikro-Inspection Tool (option).....	167
5.11	Activate QEVO Mikro-Inspection Tool (option)	168
5.12	Activating a connection to the navigation system	168
5.13	Preparing the device for the use of laser micromanipulators.....	169
5.14	Device function tests	171
6	Operation	173
6.1	Safety during operation.....	173
6.2	Software configuration.....	175
6.2.1	Touchscreen	175
6.2.2	Configure microscope	175
6.2.3	Configuring the button assignment of the hand grips, foot control panel and rocker foot switch ..	180
6.2.4	Configure stand.....	184
6.2.5	Configure displays	193
6.2.6	Configure audio	193
6.2.7	Configuring Photo	194
6.2.8	Configuring the video camera	194
6.2.9	Configuring the video recording	196
6.2.10	Configure fluorescence.....	197
6.2.11	Configure MultiVision.....	199
6.3	Manage User	201
6.3.1	Adding a user	201
6.3.2	Configuration of users	202
6.3.3	Export user settings to a USB storage medium.....	202
6.3.4	Importing the user settings from a USB storage device	203
6.3.5	Changing user	203
6.3.6	Copy a user	204
6.3.7	Delete User	204
6.4	Patient administration	205
6.4.1	Adding a patient	205
6.4.2	Changing patients	205
6.4.3	Edit patient data.....	205
6.4.4	Searching for a patient	206
6.4.5	Sorting patients.....	206
6.4.6	Export Patient Data	206
6.4.7	Importing patient data from the USB port or network	208
6.4.8	Importing patient data from the DICOM Worklist.....	209
6.4.9	Deleting patients	210
6.5	Recording video files (option)	210
6.6	Additional recording of video files on a USB storage medium or network drive (parallel recording).....	213
6.6.1	Parallel recording of video files on a USB storage medium.....	213
6.6.2	Parallel recording of video files on a network drive.....	213
6.7	Recording video files via SmartRecording.....	214
6.8	Streaming video files	216
6.9	Taking photos with the integrated video camera	217
6.10	Taking photos with the DSLR camera (option)	218

6.11	Patient data, photographs and videos	219
6.11.1	Opening the patient directory	219
6.11.2	Display of media objects (files)	219
6.11.3	Viewing patient images	219
6.11.4	Viewing patient videos	220
6.11.5	Editing patient videos	221
6.12	DICOM (option)	222
6.12.1	Importing patient data from the DICOM Worklist	223
6.12.2	Exporting patient data to DICOM	223
6.13	Use QEVO Micro-Inspection Tool (option)	224
6.13.1	Switch QEVO on / off	224
6.13.2	Display QEVO user interface	225
6.13.3	Switching the automatic switch-off function on/off	226
6.13.4	Perform white balance	226
6.13.5	Setting the image brightness	227
6.13.6	Change QEVO alignment	228
6.13.7	Show / hide marker	231
6.13.8	Switching the light on/off	231
6.14	BLUE 400 (option)	232
6.14.1	Configuring BLUE 400	233
6.14.2	Carry out BLUE 400 functional inspection	234
6.14.3	Starting the BLUE 400 fluorescence application	235
6.15	YELLOW 560 (option)	235
6.15.1	Configuring YELLOW 560	236
6.15.2	Procedure for YELLOW 560 function test	237
6.15.3	Starting the YELLOW 560 fluorescence application	238
6.16	Using the device as a hotspot	239
6.17	Connecting an external device to the device via web interface (option)	240
6.18	Using ZEISS Smart Services	241
6.18.1	Activating ZEISS Smart Services	241
6.18.2	Send Diagnostic Package	241
6.18.3	Allow Remote Service Access	241
6.19	Powering the device down	242
7	Cleaning and disinfection	243
7.1	Safety during cleaning and disinfection	243
7.2	Contamination of the system	243
7.3	Cleaning and disinfecting agent	244
7.3.1	Cleaning agents	244
7.4	Cleaning	244
7.4.1	Cleaning optical surfaces	244
7.4.2	Cleaning the monitor (touchscreen)	244
7.4.3	Cleaning mechanical surfaces	245
7.4.4	Fogging of optical surfaces	245
7.5	Disinfection	245
7.5.1	Disinfecting the mechanical surfaces	245

8	Maintenance	247
8.1	Safety during maintenance.....	247
8.2	Maintenance schedule for the operator.....	247
8.2.1	User maintenance intervals.....	247
8.3	Maintenance schedule for the authorized service.....	248
8.3.1	Service maintenance intervals.....	248
8.3.2	Performance of safety inspections	248
8.4	Remote service	249
9	Troubleshooting.....	251
9.1	Localize faults.....	251
9.1.1	React to disruptions with notifications.....	251
9.1.2	Export Log Files	251
9.1.3	Service information	251
9.2	Remote service	252
9.3	Malfunctions without notifications.....	252
9.4	Basic Function Mode following system control failure	254
9.5	Device in the de-energized state.....	255
9.6	Automatic lamp change.....	255
9.7	Manual lamp change	256
9.8	Replacing the lamp container	257
10	Technical specifications.....	261
10.1	Regulatory information	261
10.2	Radio Frequency Approval.....	261
10.3	Radio modules.....	261
10.3.1	Bluetooth module	261
10.3.2	WLAN module (option).....	262
10.3.3	RFID module.....	262
10.4	Electrical data	263
10.5	Light source.....	263
10.6	Aiming beam laser.....	264
10.7	Integrated HD camera	264
10.8	4K camera	264
10.9	Digital video outputs	264
10.10	Digital video inputs.....	265
10.11	4K video outputs	265
10.12	Audio input	265
10.13	Mechanical data	265
10.13.1	System data	265
10.13.2	Dimensions of rollers.....	266

10.13.3	Component weights.....	266
10.14	Optical data.....	267
10.14.1	Surgical microscope.....	267
10.14.2	Widefield eyepiece (magnification factor 10x)	267
10.14.3	Widefield eyepiece (magnification factor 12.5x)	268
10.15	Ambient requirements for operation	268
10.16	Ambient requirements for transport and storage	268
10.17	Operation with intraoperative MRT systems	268
10.18	Dimensional drawing.....	269
10.19	Guidelines and manufacturer's declaration for electromagnetic compatibility	270
10.19.1	EMC - Electromagnetic compatibility as per IEC 60601-1-2: 2007 (3rd Edition)	270
10.19.2	EMC - Electromagnetic compatibility as per IEC 60601-1-2: 2014 (4th Edition)	274
10.20	Security Concept.....	277
10.20.1	Objective.....	277
10.20.2	System overview.....	277
10.20.3	General principles.....	278
10.20.4	Cybersecurity relevant functions	279
10.20.5	Ensure trusted system setup	284
10.20.6	Protection of critical device functionality	285
11	Accessories and components	287
11.1	Video components.....	287
11.2	QEVO	287
11.3	Components for hybrid visualization	287
11.3.1	Main tube	287
11.3.2	Eyepieces for main tube	288
11.3.3	Co-observer tube.....	288
11.3.4	Tube for left/right co-observation	288
11.3.5	Eyepieces for left/right co-observation	288
11.3.6	Tube for face-to-face.....	288
11.3.7	Eyepieces for face-to-face.....	289
11.3.8	Tube for posterior fossa.....	289
11.3.9	Eyepieces for posterior fossa	289
11.3.10	Canon EOS DSLR Adapter Kit.....	289
11.3.11	Mouth switch.....	289
11.3.12	Magnification Changer.....	289
11.3.13	Adaptation of laser micromanipulators	290
11.4	Tube covers	290
11.5	Fluorescence targets.....	290
11.6	Foot control panel	290
11.7	FCP cable/FCP backup cable	290
11.8	Rocker foot switch	290
11.9	Sterile covers/drapes	291
11.10	Country-specific cables	291
11.11	Stereo glasses.....	291

11.12	Medical loupes	291
11.13	Lamp container of the xenon light source	292
11.14	USB stick	292
11.15	USB memory medium	292
12	Disposal	293
	Glossary	295
	Index	297

1 Notes on the instructions for use

1.1 Product name

KINEVO 900 is referred to as "Device" in these Instructions for Use.

1.2 Scope of application

The present instructions for use apply to KINEVO 900 with software release 1.8 and the following identification:

- Reference number: 6640

1.3 Purpose and storage of the documentation

These Instructions for Use explain the safety features, functions and performance parameters of the device. They contain instructions on the safe use of the device and identify measures for its care and maintenance.

Correct operation of the device is imperative for its safe and successful function.

- ▶ Read these Instructions for Use before setting up and using the device the first time.
- ▶ Keep the Instructions for Use accessible for all users at all times.
- ▶ Pass the Instructions for Use to future owners of the device.

1.4 Questions and comments

- ▶ If you have questions or comments regarding these instructions for use or about the device, contact ZEISS Service.

You can find the ZEISS contact partner for your country on the following website: www.zeiss.com/med

1.5 Conventions in this document

Certain types of information are specially marked in this document for better recognition.

1.5.1 Conventions in all text areas

- This is a list.
 - This is a second level list.

This is a cross-reference: Conventions in this document [11].

This is **highlighted text**.

This is `software code or program text`.

Names of software dialogs, fields or menus, and software messages are marked by quotation marks:

- "View" menu.
- "Do you want to save the settings?"

The steps in menu and file paths are separated by slashes:

- "File / Save as"
- "My documents / Documents"

Keys, buttons, knobs, levers and other operating controls are marked by square brackets:

- [START] key
- [Next] button

1.5.2 Conventions in a course of action

WARNING!

This is warning information about hazards that can cause death or severe injuries if not avoided.

The warning message names the possible consequences.

- ▶ This is a measure with which hazards can be prevented.

CAUTION!

This is warning information about hazards that can cause injuries if not avoided.

The warning message names the possible consequences.

- ▶ This is a measure with which hazards can be prevented.

NOTE

This is warning information about hazards that can cause property damages if not avoided.

The warning message names the possible consequences.

- ▶ This is a measure with which hazards can be prevented.

- ☑ This is a requirement that must be met before the start of a sequence of actions.

1. This is a command.

2. **CAUTION! This is a warning message about hazards that can occur during a single action.** This is a command.

⇒ This is the result of a sequence of actions.

1.6 Other applicable documents

Document type	Document title	Document number
Training record	Briefing report for newly installed systems	G-30-1714
User manual	INFRARED 800 with FLOW 800 Option	G-30-1956
User manual	Mouth switch (option)	G-30-1469
User manual	14-function foot control panel, wired (FCP) or wireless (FCP WL) (option)	G-30-1706
User manual	QEVO and QEVO ECU (option)	G-30-1976
Reprocessing instructions	QEVO	G-30-1989
User manual	Loupe EyeMag Pro S	G-30-1688
User manual	SMARTDRAPE	G-30-1996
Product overview	ZEISS Video Accessories	G-30-1888
Conformance statement	DICOM Declaration of Conformity	G-30-1952
License agreements	MICROSOFT SOFTWARE LICENSE TERMS for WINDOWS 10 IOT ENTERPRISE & MOBILE	accessible via the device's user software
Release Notes	Release Notes for ZEISS KINEVO 900	RN-30-9063

Empty page, for your notes

2 Safety notes

2.1 Target group

CAUTION!

Only to be operated by trained personnel!

These Instructions for Use are intended for physicians, nurses and other medical and technical staff who prepare, operate or maintain the device after appropriate training. It is the duty of the device owner/operator to train and brief all the operating personnel.

- ▶ Initial instruction in preparation, operation, warnings/hazards, emergency operation and transport of the device shall be provided according to Training record G-30-1714 in connection with the present Instructions for Use.
- ▶ Any further training and instruction of planned operating personnel shall be performed by the operator of the device based on these Instructions for Use.

2.2 Area of use

2.2.1 Intended use

The KINEVO 900 is a surgical microscope intended for the illumination and magnification of the surgical area and for the support of visualization in surgical procedures.

CAUTION!

Injury to the patient's eye!

The device must not be used for ophthalmological procedures.

- ▶ Make sure that no xenon light and no laser radiation enters the patient's eyes.

Optional functionalities

The BLUE 400 option is designed for a surgical microscope and can be used in tumor surgery to visualize tumor tissue marked with fluorescence dyes.

The YELLOW 560 option is an accessory for the intraoperative illumination of the surgical area with light in the wavelength range from 460 nm to 500 nm. Light in the wavelength range between 540 and 690 nm is particularly highlighted.

2.2.2 Indications

Information on the indication of use, contraindications and the target patient group is set out in the following table.

Indication of use	
Discipline:	■ Neurosurgery and spine surgery ■ Orthopaedic surgery ■ Oral and maxillofacial micro-surgery ■ P&R procedure ■ Not applicable for ophthalmology
Indication:	The ZEISS KINEVO 900 is a surgical microscope intended to illuminate and magnify the surgical area and to support visualization in surgical procedures.
Contraindication:	Provided that ZEISS KINEVO 900 is used within its indication of use, there are no known contraindications.
Patient target group:	No restriction
Age:	No restriction
Gender:	All

2.2.3 Normal use

The complete system comprises a surgical microscope with a highly flexible floor stand containing the electronics, the light source, a graphic touchscreen with video display and an optional second monitor (2D or 3D). The stand supports the operation of the system in the following combinations:

- as an optical system
- as a digital system
- as a combined optical-digital ("hybrid"), robotic visualization system

The visualization system can use mono or stereo video functionalities combined with or without optical observation (tubes/eyepieces). ZEISS KINEVO 900 is ideally suited for cranial and spinal applications in neurosurgery as well as ENT applications in the area of the auditory nerve and the base of the skull. Further fields of application include P&R procedures and oral and maxillo-facial surgery. The system is also ideally suited to multidisciplinary use in microsurgery.

It is also designed for surgical interventions in which an endoscope

and a surgical microscope are used at the same time, particularly with its microvisualization tool which can be added as an optional extra and simply plugged into the system so that the surgeon can detect hidden details or "see around the corner".

The system is equipped for the connection of navigation systems, intraoperative neuromonitoring systems and for the data communication with external network systems. Optional fluorescence modules enable a better visualization of tumor and healthy tissue as well as better vascular visualization.

The system is intended for use in hospitals, clinics or other medical facilities.

The functions of the surgical microscope (OPMI) and the stand are controlled via the central control unit in the console. The interactive graphic touchscreen enables the user to configure all necessary settings. The user may trigger these functions using the buttons on the handgrips or on an optional foot control panel.

The system may be operated only by physicians, nurses and other medically trained OR staff who have received appropriate training and observe the instructions for use.

The installation conditions and the use of the system must meet microsurgical requirements:

- Low vibration
- Dust-free environment
- Level, horizontal positioning
- Avoidance of extreme mechanical stress
- EMC-specific environment based on a professional healthcare facility (use of RF surgical devices in combination with the microscope)

2.2.4 Possible damage/injury from burns due to the high light intensity

CAUTION!

Tissue damage and burns at high light intensity!

A high light intensity over a longer period of time can damage tissue and cause burns in the area of the illuminated surgical field.

- ▶ Read and observe the factors which can increase the risk of burns described below.
- ▶ Read and observe the recommendations on how you can reduce the risk of burns.
- ▶ Make sure that the "Focus Light Link" function is always switched on.

⇒ The "Focus Light Link" function limits the maximum adjustable light intensity in the focal plane depending on the working distance. Possible tissue damage due to an unintentional high light intensity can thus be avoided. This function should only be deactivated if a higher light intensity is required for the current application.

General

- The device is equipped with a powerful xenon illumination.
- If used improperly, excessive illumination intensities may lead to third-degree burns.
- A device with light source switched on must never be left unattended.

Various factors contribute to the risk of burn injuries:

System-related factors

- If a high magnification is used, the diameter of the field of view and the light intensity at the surgeon's eye and on the integrated cameras both decrease whereas the light intensity in the surgical field remains the same. This effect is increased by the use of certain components, such as e.g. eyepieces with higher magnification or the foldable tube f170/f260 with tube magnification (PROMAG function). When working at maximum magnification, you should therefore pay particular attention to the set light intensity to prevent burns, especially of the surrounding tissue.

Surgery-related factors

- The size of the luminous field influences the risk of injury in two respects: With a large luminous field diameter, areas of the skin are illuminated which are not monitored as strictly by the surgeon and which are not moistened sufficiently. These areas are particularly at risk of injury. These injuries can be prevented by adjusting the illuminated-field diameter to the smallest size required for the relevant operation.
- If the luminous field is reduced in size, the intensity increases because the light is focused more. So, if possible, the intensity should be lowered as soon as the size of the illuminated field is reduced. The device can optionally be equipped with Automatic Light Field Limitation.
- A long surgical procedure increases the risk of injury, in particular if a standard procedure takes considerably longer than usual.
- Injuries in the peripheral area can be prevented by covering this area with moist, sterile gauze. Moisten the gauze at regular intervals to prevent the area from drying out or heating up. The risk is increased if dry drapes are used to cover such areas.
- Some areas of the patient's body may be more sensitive than others.
- Certain preparations of the surgical field, local vasoconstrictive medications and drapes at the wound edge may also result in a higher risk of injury (drapes may heat up to varying degrees depending on their color and moisture content).
- When working in deep channels, the working distance from the illuminated wound edges is shorter than the distance from the actual surgical field. If you increase the light intensity to better illuminate a deeper surgical field, the risk of burns increases for the wound edges.

Patient-related factors

- The general condition of a patient's health may contribute to the risk of injury.
- The skin type may be decisive for the risk of burns or injuries caused by light.
- Certain medications (e.g. light sensitizing substances) also affect the sensitivity to light.
- The interaction of heat and antimicrobial substances in incision foils may lead to an increased reaction of the patient to these substances.

Recommendations

- Make sure that the "Focus Light Link" function is switched on. The "Focus Light Link" function limits the maximum adjustable light intensity in the focal plane depending on the working distance. Possible tissue damage due to an unintentional high light intensity can thus be avoided. This function should be deactivated temporarily at most and only if a higher light intensity is required for the current application. Therefore, please note the recommendations for preventing burns.
- The initial illumination intensity should be preset to a low value.
- If used, drapes should also be remoistened at regular intervals in order to prevent heat from accumulating underneath the drape.
- The risk of burns can be reduced by constantly irrigating the illuminated surgical field and by keeping it moist.
- Using the buttons on the hand grip or foot control panel, the surgeon can set the illumination intensity to the value required for the procedure. Please note that the intensity increases with decreasing luminous field size if the spot function is used. For this reason, the intensity should be set only after the size of the illuminated field has been changed.
- The device features an "Automatic Light Field Limitation" which should not be turned off.
- Reduce the luminous field diameter when working in deep channels, so that the wound edges are not illuminated at a high light intensity wherever possible.
- If the magnification factor is increased during surgery, the light intensity on the observer's side decreases. For balancing, the illumination intensity must be increased. This loss in image brightness is automatically compensated if the "Auto Brightness" function is activated.
- Never leave a device unattended when its light source is switched on.
- Switch off the light when the microscope is not being used, and ensure it is not directed at unprotected bare skin.

2.2.5 Electromagnetic compatibility

The device is subject to specific requirements with regard to electromagnetic compatibility (EMC).

Take the following precautionary measures to avoid EMC malfunctions:

Procedure

- ▶ Observe the Instructions for Use.
- ▶ Observe the EMC guidelines in the section "Technical specifications".
- ▶ Only use accessories, cables and spare parts which have been approved by ZEISS for this device.
- ▶ If you use radio equipment or wireless transmission components: Maintain a minimum distance of 30 cm from all components of the device.
- ▶ If you install the device in the vicinity of other devices or stack it with other devices: Check whether normal operation is possible in this arrangement.

2.3 Responsibilities and duties of the operator

Operating personnel

The device may only be operated by instructed and trained personnel.

- ▶ Ensure that the operating personnel have been trained and instructed.
- ▶ Ensure that the operating personnel have read and understood the instructions for use.
- ▶ Keep the Instructions for Use available at all times for the operating personnel.
- ▶ To simplify access for all operating personnel, order additional copies of the Instructions for Use from ZEISS as required.
- ▶ Define the required skills for handling the device and provide information on who is authorized for which activities.
- ▶ Define rules for reporting errors and damage, and make sure that people are informed of these. Notification to manufacturers and authorities [25]
- ▶ Provide the necessary protective clothing.
- ▶ Regularly check compliance with the national laws and regulations concerning accident prevention and occupational health.

Safety inspections

- ▶ To prevent a decrease in device safety due to aging and wear: Have a safety inspection in compliance with IEC 62353 performed on time and to the prescribed extent. Observe the applicable national regulations.

The safety inspections may only be performed by the manufacturer or qualified personnel.

- ▶ Comply with the specified time limits.
- ▶ Perform the inspections to the prescribed extent.

The device should, as a minimum, undergo the following checks and safety inspections:

- Check that the Instructions for Use are present
- Visually inspect the device and its accessories for damage, as well as legibility of markings and labels
- Leakage current test
- Test of protective ground conductor
- Function and wear test of the brakes
- Function test of all switches, buttons, sockets and LEDs on the device

Maintenance and inspection

- ▶ To ensure safe operation of the device and reach the expected service life: Observe the maintenance and inspection intervals that are specified in these instructions for use.

Changes to the product

- ▶ **WARNING:** This device must not be modified without the manufacturer's permission. If the device is modified, suitable inspections and tests must be performed to ensure that the device can still be used safely.

Accessories and additional equipment

- ▶ If you wish to connect accessories or additional equipment to the device: Contact your ZEISS representative [11].

Additional equipment connected to medical electrical devices must demonstrably comply with the corresponding IEC or ISO standards (e.g. IEC 60950 for data processing equipment).

Furthermore, any configurations must comply with the requirements stipulated in standards for medical systems (see IEC 60601-1).

If you connect additional devices to medical electrical systems, you are considered to be a system configurer and are thus responsible for ensuring that the system complies with the normative requirements for systems.

Local laws take precedence over the above normative requirements.

Dangers arising from connection to a network

Before connecting the device to a network, observe the following safety measures to prevent injury or damage:

- ▶ The user (or IT manager) is responsible for ensuring that no computer viruses are transmitted to the device via the network connection.
- ▶ Ensure a data transmission rate of 1 Gbit/s (Fast Ethernet) and the conformity of your network configuration with internet protocol IPv4 so that patient data can be exported to your network reliably, safely and error-free.

Changes to the network

The following subsequent changes to the network may result in new risks:

- Changes to the network configuration
- Connecting additional devices to the network
- Disconnecting devices from the network
- Updates of devices connected to the network
- Upgrades of devices connected to the network
- ▶ Analyze and eliminate any new risk factors arising through such changes to the network.

2.3.1 Requirements for the network connection

The following requirements apply to the network in which the ZEISS KINEVO 900 is integrated:

LAN

Name	Value
Ethernet protocol IPv4	TCP/IP
Video Streaming	RTP via UDP
Gigabit Ethernet	1000BASE-T, IEEE 802.3 Paragraph 40
Assignment of IP addresses	Static IP or DHCP (configurable within the device)
LAN connection to the navigation system	
Latency	< 3 ms
Data transmission rate	> 1 Mbit/s
Parallel video recording for release in the network	
Data transmission rate (without SmartRecording)	> 40 Mbit/s

Name	Value
Data transmission rate (with SmartRecording)	> 80 Mbit/s
Live streaming of camera image	
Data transmission rate	> 40 Mbit/s

WLAN

Name	Value
Safety standard	WPA2/PSK or WPA2/PEAP
WLAN standard	802.11a/b/g/n
Ethernet protocol IPv4	TCP/IP
Video Streaming	RTP via UDP
Assignment of IP addresses	Static IP or DHCP (configurable within the device)
Parallel video recording for release in the network	
Data transmission rate (without SmartRecording)	> 40 Mbit/s
Data transmission rate (with SmartRecording)	> 80 Mbit/s
Live streaming of camera image	
Data transmission rate	> 40 Mbit/s

WLAN-Hotspot

Name	Value
Safety standard	WPA2/PSK or WPA2/PEAP
WLAN standard	802.11a/b/g/n
Ethernet protocol IPv4	TCP/IP
Video Streaming	RTP via UDP
Assignment of IP addresses	DHCP
Live streaming of camera image	
Data transmission rate	> 40 Mbit/s

Network configuration

Name	Value
Assignment of IP addresses	Static IP or DHCP (configurable within the device)
Shared network drive	SMB protocol (Server Message Block), also known as CIFS (Common Internet File System)
Streaming client (via LAN, WLAN, WLAN hotspot)	
Log	RTSP (RTP via UDP)
2D video codec	■ Encoding: H.264/MPEG-4 AVC / High Profile / Level 4.2 ■ HD resolution: 1920x1080p50/60 ■ LowRes resolution: 960x540p25/30
3D video codec	■ Encoding: H.264/MPEG-4 AVC / High Profile / Level 4.2 ■ Resolution: 1920x1080p50/60 TopBottom

2.3.2 Messages to manufacturer and authorities

If a serious incident occurs in connection with this medical device affecting the operator or another person, the operator (or person responsible) must report this serious incident to the manufacturer of the medical product. In the European Union, the operator must report this serious incident to the competent authority in his/her country.

2.4 Measures and duties of the operator

Electrical safety

- ▶ Always switch off the device before connecting it to or disconnecting it from the power supply, for cleaning its surface, or if it will not be used for a prolonged period of time.
- ▶ Only connect the device to a power supply that complies with the values specified on the rating label.
- ▶ Do not use multiple sockets!
- ▶ Do not use extension cables!
- ▶ Do not touch the device if your body is electrostatically charged and the device is not grounded.
- ▶ Connect the device via the potential equalization connection (according to IEC 60601-1) to other active devices with the same ground potential or connect it to a protective ground connection.
- ▶ Please observe the information on electromagnetic compatibility (EMC).

The device contains freely accessible live components. If you remove the housing, you run the risk of electric shock.

- ▶ Never open the device!

Environmental conditions

- ▶ Make sure that the installation conditions and the operation of the device comply with the surgical requirements:
 - Low vibration
 - Clean environment
 - Avoid extreme mechanical stress
- ▶ Do not use power-operated devices included in the delivery package
 - in explosive atmospheres,
 - at a distance of less than 25 cm from flammable anesthetics or volatile solvents such as alcohol, benzine or similar substances.
- ▶ Do not use or store the device in damp rooms. Do not expose the device to water splashes, dripping water or sprayed water.
- ▶ Ensure that fluids cannot enter the device.

Symbols and labels

- ▶ Note the symbols and labels attached to the device!

Transport

- ▶ Only transport the device over long distances (e.g. relocation, return for repair) in its original packaging or special return packaging.
- ▶ Please contact your dealer or ZEISS Service for this purpose.

2.5 Liability and warranty

The warranty and liability depend on the contractually specified conditions.

Do not modify the device without permission.

- This device must not be modified without the manufacturer's approval. If the device is modified, suitable inspections and testing must be completed to ensure that it can still be used safely.
- The manufacturer is not liable for damage caused by unauthorized use of the device. Furthermore, this will forfeit any rights to claim under warranty.

3 Device description

3.1 General

ZEISS KINEVO 900 is a robotic visualization device for demanding neurosurgical applications. A modular hard and software concept enables the successive expansion of the system with additional options such as co-observation lenses, fluorescence modules or the preparation for the connection of navigation systems.

The fully digital configuration enables exclusively digital operation with an external 3D video monitor. Optical co-observation cannot be used in this configuration.

Only options that have been installed and activated with a corresponding license are displayed on the touchscreen.

3.2 Innovative functionalities

3.2.1 Kinematic and robot-assisted positioning functions

KINEVO 900 offers the user new and innovative movement and positioning functions for simplified and more precise pre- and intra-operative microscope positioning.

- Manual "PointLock" movement mode

The [SB] (Selected Brakes) buttons can be assigned using the PointLock function. This function enables the microscope to be moved around a constant visual field centre and repositioned. The focus is at the same time automatically adjusted from 200 to 625 mm within the available focal length range.

- Motorized "PointLock" XY movement mode

This new motorized XY movement function enables the adjustment of the microscope while maintaining the image centre and the working distance. The motor control can be triggered via the XY joystick on the handgrip or hands-free with the FCP.

- Position Memory

The current position of the microscope as well as the working distance and the magnification can be saved at any time. Via the handgrip and the FCP or the user interface on the touchscreen, these stored positions can be called up again and moved to with high accuracy. Important intraoperative landmarks or anatomical details can be called up again and set in this way.

- **Motorized "Microscope" XY movement mode**

The microscope can be moved by motor in the three axes of its suspension, 4, 5 and 6. This enables motorized tilting and swiveling of the focal point in the XY direction with a fixed working distance. The stand does not move in the process.

- **Motorized "Stand" XY movement mode**

The microscope can be moved by motor in the three axes 1, 2 and 3 in the XY focal plane without tilting or swiveling; the alignment of the eyepieces always remains the same (e.g. horizontal). This enables precise motorized XY movement of the focal point at a fixed working distance. Only the axes of the stand move in the process.

- **Motorized movement to the system park position**

The device can be automatically moved to its park position.

- **Motorized movement to a drape position**

A drape position is saved as a factory setting, however, can be changed arbitrarily. This drape position can be redefined and saved by the user in order to enable more efficient draping always performed in the same position.

3.2.2 QEVO handheld Micro-Inspection Tool (option)

QEVO (option, see Instructions for Use G-30-1776) is a new type of hand-held, digital and completely sterilizable exploration tool for the display of anatomic details not accessible with the microscope (e.g. a view behind an aneurysm to be clipped). It can be plugged into the console of the KINEVO 900 easily and quickly as required, thus replacing an additional endoscope system which would require both time and effort to prepare. The clear high-resolution HD video image of QEVO is displayed on the system monitor or an externally connected video monitor by pressing a button on the handgrip or the FCP.

3.2.3 Integrated 3D video system for observation without eyepieces (option)

ZEISS KINEVO 900 can be used with a fully integrated stereo video camera as a purely digital visualization system. The observation is carried out in this case using stereo glasses via an external 3D video monitor or an additional 3D video monitor (component of the stereo option), which is installed directly on the console instead of the second HD system monitor. If configured accordingly - depending on the surgical situation at hand and individual preference - the system can also be used with a tube and eyepieces and/or stereo glasses and monitor.

NOTE

Application of BLUE 400 only with eyepiece-based optical observation!

Tubes and eyepieces must be used when the BLUE 400 option is applied.

3.2.4 Data management

- **Data storage in a shared network (option)**

The KINEVO 900 can be connected to the clinic network via a LAN connection (and optionally via WLAN, see below). A correspondingly configured server directory enables fast and easy data transfer of images and videos to this network drive for central data storage and fast accessing of this data from every PC/laptop connected to the network for editing and processing (e.g. with a video editor). This option also enables parallel video recording (option) on this network drive and prevents e.g. time-consuming data transfer using USB media.

- **Wireless WLAN data transfer (option)**

An optional WLAN module with a WiFi hotspot function simplifies the management of images and videos stored with the KINEVO 900 enormously. With this option, a wireless connection to the clinic network or server directory can be configured instead of a LAN connection for time-saving data transfer. A WiFi hotspot for connections to mobile devices (e.g. tablets or smartphones) can be additionally activated and used, thus enabling users to access the patient data directory in the KINEVO 900 and download the desired images and video clips quickly and easily. All network, WLAN and WiFi connections are password-protected.

3.2.5 Navigation functionality

The KINEVO 900 is prepared for the display of navigation information and the connection of external navigation systems (factory-integrated standard navigation interface and MultiVision system). In connection with an extended navigation license, the new robotic positioning capabilities featuring motors in all system axes permit automatic alignment of the microscope to a trajectory predefined in the navigation planning software, provided that the connected navigation system supports this function.

3.2.6 Video recording functionality (option)

The KINEVO 900 offers the user a considerably more effective video recording functionality. When using a stereo video camera, the user can choose between a 2D or a 3D recording of the surgical procedure. In addition to the normal HD video recording (activated via hand grip, FCP or touchscreen user interface), a disk space optimized recording with a reduced resolution also can be activated, e.g. in order to record and save daily surgical procedures in their full length. An intuitive video editor with a cut & merge function enables fast and easy video editing on the system prior to downloading. A video streaming option can be activated if the system is connected to the clinic network. This option makes it possible to transmit the surgical procedure via LAN/WLAN e.g. to a remote auditorium with any accessible IP address or to a connected mobile device via WiFi.

Note:

If mobile systems are to be used for the display of the video stream, ZEISS recommends using the VLC media player (VideoLAN) for Apple IOS and Windows.

3.3 Device marking

CAUTION!

Risks resulting from non-observance of symbols and markings

Failure to observe warning symbols and markings on the product may result in misuse of the device, may damage the device, and may endanger the health of the patient and the operator.

- ▶ Before operating the device, familiarize yourself thoroughly with the symbols and markings on the product.
- ▶ Observe the symbols and markings on the device at all times!

3.3.1 Labeling on the microscope


Figure 1: Labeling on the microscope

Pos.	Symbol	Explanation
1		Max. load: 6 kg
2		Open, right/left co-observation port
3		Warning of laser beam (labeling for autofocus option)
4	Varioskop	Brand name
5		Brand / ZEISS logo
6	Z	Zoom
7	F	Focus
8	302584-9000-888 XXXXXXXXXX	Microscope identification label (inside)
9		Switch over, right/left co-observation port
10		Luminous-field diaphragm
11		Focus/working distance
12		Adjust zoom magnification manually

3.3.2 Approval label / UDI label

Figure 2: Approval label / UDI label

Pos.	Explanation
1	Approval label ■ Manufacturer: Carl Zeiss Meditec AG ■ Manufacturer symbol: ■ Manufacturer's address ■ Device name: KINEVO 900 ■ Reference and serial number ■ Date of manufacture symbol: ■ Date of manufacture, year-month-day ■ Rated voltage: 100V - 240V ■ Max. connected load: 1350VA ■ Line frequency range: 50Hz - 60 Hz ■ Designation of origin for industrial products: MADE IN GERMANY ■ WEEE mark: ■ CSA mark: ■ CE label ■ Maximum total mass: Max. 395 kg
2	UDI label ■ Machine-readable label (barcode) ■ UDI Device Identifier (UDI-DI) ■ UDI Production Identifier (UDI-PI)

3.3.3 Labeling on stand, Part 1

Figure 3: Labeling on stand

Item	Symbol	Explanation
1	KINEVO 900	Device name
2		Brand / Logo
3		Seal for drape suction
4		Risk of crushing!
5		Transport position: ■ Observe the Instructions for Use ■ Pushing position ■ General warning

3.3.4 Labeling on stand, Part 2


Figure 4: Markings on the stand, part 2

1		Standby/ON-OFF switch
2		RFID scanner or reading device
3	USB OPEN	Open USB cover Note: This marking is only present on older devices [46].
4		4K video output ■ Number 1: 4K HDMI output ■ Number 2: 4K HDMI output
5		Video input ■ Number 3: Display port input ■ Number 4: Display port input ■ Number 5: HDMI/DVI input

Pos.	Symbol	Explanation
6		Video output ■ Number 6: HDMI/DVI output ■ Number 7: Display port output ■ Number 8: Display port output ■ Number 9: HD-SDI / 3G-SDI output ■ Number 10: HD-SDI / 3G-SDI output
7		Remote port
8		Foot control panel, 14 functions
9		Audio in
10		Audio out
11	LAN	Ethernet port (LAN)
12		Rocker switch
13		4K video output (present only on 4K option, not yet available) ■ Number 11: 3G-SDI output, left (2D/3D) ■ Number 12: 3G-SDI output, left (2D/3D) ■ Number 13: 3G-SDI output, left (2D) ■ Number 14: 3G-SDI output, left (2D) ■ Number 15: 3G-SDI output, right (3D) ■ Number 16: 3G-SDI output, right (3D)

3.3.5 Labeling on stand, Part 3




Figure 5: Labeling on stand, Part 3

Pos.	Symbol	Explanation
1		Risk of crushing!
2		Approval label for INFRARED 800 (option) provides information on: ■ Manufacturer ■ License data ■ UDI labeling ■ CE label
3		UDI label for ZEISS KINEVO 900, provides information on: ■ UDI labeling ■ Serial No. ■ Device name
4		Fuse display F1

Pos.	Symbol	Explanation
5		Read and observe the Instructions for Use.
6		Non-ionizing electromagnetic radiation
7		Potential equalization
8		Approval label for ZEISS KINEVO 900 provides information on: ■ Manufacturer ■ System data ■ Power data ■ UDI labeling ■ Type of disposal ■ total weight ■ CSA label ■ CE label
9		SIP label, provides information on: ■ Manufacturer ■ Manufacturer's contact information ■ SIP number of device (SIP = Service Identification Program)
10	6640xxxxxx	Serial number of device
11	MD	MD labeling (Medical Device)* * Depending on the respective country-specific registration status
12	NAV	Port for navigation system
13		Identification label, options

Pos.	Symbol	Explanation
		Identification for installed options

3.3.6 Labeling on stand base


Figure 6:

Item	Symbol	Explanation
1		Pedal locking tabs
2		Straight-ahead travel pedal

3.3.7 Labeling on light source



Figure 7: Labeling on light source



Pos.	Symbol	Explanation
1		Lamp replacement/manual lamp replacement (removal of light source cover)
2		Risk of crushing!
3		Hot surface
4		General Warning
5		"Observe disposal regulations" label Do not dispose of electrical or electronic devices along with normal domestic waste.
6		Manual lamp change
7		Identification label, lamp 300 W+ xenon
8		Light 1

Pos.	Symbol	Explanation
9		Light 2
10		Identification label, "HLQ 300 Future" light source

3.3.8 Labeling for autofocus option

Figure 8: Labeling for autofocus option

Pos.	Symbol	Explanation
1		Warning of laser beam
2		Country-specific labels: Labeling in accordance with IEC 60825-1:2007 and IEC 60825-1:2014 Laser radiation - Do not look into the beam or view directly with optical instruments Performance features ■ Laser class: II ■ Power output max: < 1 mW ■ Wavelength: λ 635 - 645
3		USA label: Laser class 2

3.3.9 Labeling on packaging

Symbol	Symbol	Explanation
	Indication of direction "This side up"	Indicates the correct upright position of the package.
	Fragile	Handle with care
	Keep dry	Protect packaging and packaged contents from wetness.
	Do not stack	Stacking of the packages is not permitted. No load should be placed on the package.
	Permissible temperature	The product may only be transported and stored at a temperature range of min. -20°C to max. +60°C.
	Packing unit	Specification of the number of packing units
	Permissible relative humidity	The product may only be transported and stored at an humidity of min. 10% and max. 90% RH.
	Permissible atmospheric pressure	The product may only be transported and stored at an atmospheric pressure of min. 500 hPa and max. 1060 hPa.
KINEVO 900		Device name
		Brand / Logo
		Min. 1.2 m forklift length
	Center of gravity	Indicates the package's center of gravity.

Symbol	Symbol	Explanation
		Shockwatch
		Tiltwatch

Table 1: Symbols for transport and storage

3.3.10 Displaying radio frequency approvals on the monitor

Action

1. Tap on  Settings → **i** Info.

2. Scroll down in the "Info menu": Open → Radio Frequency Approval.
⇒ The markings of the existing radio frequency approval are displayed.
3. Scroll down in the "Radio Frequency Approval" display.
4. Close the "Radio Frequency Approval" display by tapping on the black surface next to the display.

3.3.11 Display Windows EULA on the monitor

Action

1. Tap on  Settings →  Info.

2. Scroll down to the section titled "Windows End-user license agreement (EULA)"
3. Tap on the button [Open].
⇒ The Windows End-user license agreement is displayed.
4. Close the display by tapping on the black surface next to the display.

3.4 Structure of the device

3.4.1 Device overview


Figure 9: Overview of device

1	Main monitor (touchscreen):	2	Vertical arm
3	Horizontal arm	4	Microscope suspension
5	Microscope	6	Console
7	Lamp housing	8	Stand base
9	Connector panel	10	Transport handle, 2x
11	Second monitor (option)		

3.4.2 Overview of connector panel


Figure 10: Overview of connector panel

1	⏻	Powering the system up/down
2		Connector panel for external storage devices [46]
3	Ⓜ➡	Remote port (AUX) for external signal, max. 24 V / 0.5 A, galvanically decoupled as per IEC 60601-1
4	↙	Port 14-function foot control panel
5	🔊➡	Audio in
6	LAN	Ethernet port (LAN), galvanically decoupled as per IEC 60601-1
7	↙	Rocker foot switch port
8	🔊➡	Audio out
9	🗺️	Port for navigation system

Pos.	Symbol	Name
10		Potential equalization: For connecting the system to the potential equalization system in compliance with IEC 60601-1.
11		Power input socket
12		Automatic circuit breaker
13		Connection socket for hand-held, digital exploration tool QEVO

Connector panel for external storage devices

The connector panel for external storage devices varies depending on the generation of the ZEISS KINEVO 900: Depending on the build, the connector panel may have a cover which can be opened by means of a push button.


Figure 11: Connector panel for external storage devices (depending on the build)

Item	Symbol	Name
1	USB OPEN	Push button for opening the USB cover Note: The push button is only present on older devices.
2	USB	2x USB port (3.0)
3		Shelf for external USB mini HDD
4		USB cover Note: The cover is only present on older devices.

3.4.3 Overview of video connector panel

The video ports are consecutively numbered from 1 to 10 on the device.


Figure 12: Overview of video connector panel

1	4K 	Video output	4K: HDMI video output, left Output signal: HDMI 2.0 3840x2160p50/60 (ODU)
2	4K 	Video output	4K: HDMI video output, right Output signal: HDMI 2.0 3840x2160p50/60 (ODU)
3		Video input: Display port	Connection of navigation video source for MultiVision overlay
4		Video input: Display port	
5		Video input: HDMI/DVI	Video input for external video source Input signal: HDMI 1920x1080p50/60, 1280x720p50/60 (Lemo)
6		Video output: DVI	■ 2D/3D for stereo option ■ Live full screen/Endo/PIP ■ Touchscreen imaging

Pos.	Symbol	Type	Use
			Output signal: DVI-D 1920x1080p50/60 (Lemo)
7		Video output: Display port (stereo)	Camera signal video output, e.g. also for external monitor (option) Output signal: 1920x1080p50/60 (DP)
8		Video output: Display port (mono)	2D live image touchscreen, mono left Output signal: 1920x1080p50/60 (DP)
9		Video output: HD-SDI / 3G-SDI	Live camera image (corresponds to eyepiece view) with MultiVision overlay (option) Mono camera (HD): video output, left Stereo camera (3DHD): video output, left Output signal: 3G-SDI or HD-SDI 1920x1080p50/60 ("Interlace" or "Progressive", selectable), 1xBNC
10		Video output: HD-SDI / 3G-SDI	Live camera image (corresponds to eyepiece view) with MultiVision overlay (option) HD-SDI (1080i) 3G-SDI (1080p) Stereo camera (3DHD): video output, right Output signal: 3G-SDI or HD-SDI 1920x1080p50/60 ("Interlace" or "Progressive", selectable), 1xBNC

ZEISS KINEVO 900 does not support any active DP systems on
DisplayPort outputs 7 and 8.

3.4.4 Overview of video connector panel (4K 2D/3D option)*

If your device is equipped with the 4K 2D/3D option, you can connect an external 4K monitor. The ports for the 4K 2D/3D option are located on the left side of the connector panel and consecutively numbered from 11 to 16.


Figure 13: Overview of video connector panel (4K 2D/3D option)*

11		3G-SDI video output, left (2D/3D)	QuadSDI 3840x2160 50/60p (4x8NC)
12		3G-SDI video output, left (2D/3D)	
13		3G-SDI video output, left (2D)	
14		3G-SDI video output, left (2D)	
15		3G-SDI video output, right (3D)	
16		3G-SDI video output, right (3D)	

3.4.5 Overview of microscope


Figure 14: Overview of microscope

1	Left handgrip	2	Hand grip clamping
3	Microphone	4	Port for assistant tube (option)
5	Port for navigation system antenna (option)	6	Right handgrip
7	Port for co-observer right (option)	8	Microscope cover/ port for surgeon's tube (option)
9	Port for co-observer left (option)		

3.4.5.1 Hybrid visualization configuration options (option)

The information in this paragraph only applies to the optional configuration of the microscope as a system for hybrid (optical-digital) visualization.


Figure 15: Configuration options

1	Straight binocular tube	2	Tilttable binocular tube
3	Folding binocular tube	4	Tubes for face-to-face viewing (items 1, 2, 3)
5	Rotary adapter	6	Stereo co-observation module with tubes (items 1, 2, 3)
7	Adapter plate for connecting laser micromanipulator	8	Angle optics*
9	Mouth Switch	10	Photo adapter for DSLR
11	Magnification changer, 3-position*		

* When assembling the microscope, be careful not to mount the 3-position magnification changer (item 11) and the angle optics (item 8) on the microscope body at the same time. In this case, the microscope body would be too heavy for the autobalance function.

3.4.5.2 Overview of binocular tubes

The information in this paragraph only applies to the optional configuration of the microscope as a system for hybrid (optical-digital) visualization.


Figure 16: binocular tubes

1	180° tiltable tube, focal length $f = 170$ mm	2	Straight tube, focal length $f = 170$ mm
3	Foldable tube, focal length $f = 170/260$ mm	4	Foldable tube for mouth switch, focal length $f = 170/260$ mm

3.4.5.3 Widefield eyepiece (10x or 12.5x)

The information in this paragraph only applies to the optional configuration of the microscope as a system for hybrid (optical-digital) visualization.


Figure 17: Widefield eyepiece (10x or 12.5x)

1	Magnetic coupling	2	Diopter scale
3	Eyecup		

3.4.5.4 Stereo co-observation module

The information in this paragraph only applies to the optional configuration of the microscope as a system for hybrid (optical-digital) visualization.


Figure 18: Stereo co-observation module

1	Stereo coobservation module with 2 pivot joints
2	Locking lever

3.4.5.5 Digiskop cover


Figure 19: Digiskop cover

1	Cover for main observer tube	
---	------------------------------	--

3.4.5.6 Additional Light

Schematic of shadow adjustment via additional illumination

The additional illumination can be activated via the touchscreen.


Figure 20: Additional Light

1	Main illumination	2	Shaded area of main illumination
3	Surgical field	4	Special illumination to brighten shadows
5	Objective lens of surgical microscope		

3.5 Control elements and displays

3.5.1 Screen keyboard

Figure 21: Screen keyboard key assignment

Item	Name	Explanation
1	Shift key	Switches between upper-case and lower-case letters.
2	Space key	-
3	Delete key	Deletes the character to the left of the cursor.
4	Close	Ends the input and closes the screen keyboard.

3.5.2 Display: Editing

A rotating editing symbol appears during editing. A load symbol is displayed while larger amounts of data are being saved, imported or exported.


Figure 22: Progress display

During longer editing periods, the percentage display flashes inside the circle symbol.

3.5.3 Online help


Figure 23: Call-up and display of online help

Pos.	Symbol	Name	Explanation
1		Button [Help videos]	Tapping on the button opens a separate window with Help videos on various topics. Help videos feature English-language voiceover only.
2	?	Button [Online help]	Tapping on the ? button opens the Help window on the right side and the button [Help videos].
3		Help window	If the entire help text is not displayed, scroll down on the touchscreen.

The content of the Help window adapts to the respective opened menu. To close the Help window, tap on the [Online help] button again.

3.5.4 Error messages on the monitor

Information on a problem occurring during operation is detected automatically and displayed as an error message. The device saves this error message in a log file which can be exported (Settings - Extras - Export Log Files).

NOTE

Error messages on the monitor!

If a malfunction occurs during operation, a corresponding error message is displayed on the monitor.

- If an error occurs that you cannot correct with the aid of the chapter "What to do in the event of malfunctions", label the device as out of order and contact ZEISS Service.

Figure 24: Error message on the monitor

Prerequisite

- ☒ Error message appears on the monitor

Action

1. Carefully read through the error message and follow the instructions displayed therein.
2. Acknowledge the error message either directly on the monitor or by pressing the center joystick button on the **left** handgrip.
3. If the error message appears on the monitor again and the problem persists, export the log file and contact your responsible ZEISS Service organization.

3.5.5 Controls on the microscope


Figure 25: Controls on the microscope

1		Fastening screw	Attach the assembly: Tighten the fastening screw hand-tight.
2		Adjust pivoting mirror	Image outputs: Face-to-face image outputs: Right and left.
3		Fastening screw	Attach the assembly: Tighten the fastening screw hand-tight.
4		Socket (under cover)	Antenna connector (external navigation system accessory).
5		Camera release socket	For shutter release of an external camera.
6		Mouth switch socket	For connecting the optional mouth switch.
7		Adjust luminous field diameter manually	Turn CCW: Larger Turn CW: Smaller

Item	Symbol	Name	Explanation
8		Adjust focus/ working distance manually	Turn CCW: Larger Turn CW: Smaller
9		Adjust zoom magnification manually	Turn CCW: Smaller Turn CW: Larger

3.5.6 Controls of binocular tubes

The information in this paragraph only applies to the optional configuration of the microscope as a system for hybrid (optical-digital) visualization.


Figure 26: Controls of binocular tubes

Pos.	Symbol	Name	Explanation
1		Eyepiece mount	Insert widefield eyepieces into eyepiece mounting as far as they will go.
2		Adjustment wheel for eye distance	Setting the eye distance (interpupillary distance)
3		Adjustment knobs for eye distance	Setting the eye distance (interpupillary distance)
4		Eye distance scale	Display of eye distance (interpupillary distance)

Pos.	Symbol	Name	Explanation
5		Locking device for rotation function	After pressing the locking button, you can turn the tube to the left or right.
6		Adjustment wheel for integrated magnification changer (PROMAG)	Changing to a higher magnification.

3.5.7 Controls for widefield eyepieces

The information in this paragraph only applies to the optional configuration of the microscope as a system for hybrid (optical-digital) visualization.


Figure 27: Control elements on the widefield eyepiece

Pos.	Name	Explanation
1	Diopter scale	Adjustable from -8 dpt. to +5 dpt.
2	Diopter setting ring	For adjusting the diopter setting.
3	Adjustable eyecup	Can be unscrewed to prevent the entry of scattered light.

3.5.8 Controls of handgrips

Both hand grips are of the same design. You can configure the left and the right hand grip with the same functions or with different functions.



Figure 28: Button assignment of hand grips

Pos.	Name	Explanation
1	Button A , configurable	Factory setting: trigger photo
2	Button B , configurable	Factory setting: autofocus
3	Brake release button (SB)	Brake release button for selected rotary axes and PointLock (factory setting: PointLock)
4	Brake release button (AB)	Brake release button for all stand and microscope axes
5	Button D , configurable	Factory setting: increase illumination intensity Depending on the current operating mode/activated option, the yellow buttons of the right hand grip may have certain special functions (see separate description in the corresponding chapters).
6	Button C , configurable	Factory setting: reduce illumination intensity Depending on the current operating mode/activated option, the yellow buttons of the right hand grip may have certain special functions (see separate description in the corresponding chapters).
7	Rocker switch G / H	Zoom + / Zoom - (factory setting)

Pos.	Name	Explanation
8	Rocker switch E / F	Focus + / Focus - (factory setting)
9	Joystick	Motorized XY fine adjustment Depending on the current operating mode, the joystick of the right hand grip may have certain special functions (see separate description in the corresponding chapters). Messages displayed on the touchscreen can be acknowledged by pressing the center of the joystick button on the left hand grip.

3.5.9 14-function foot control panel

Pos.	Name	Explanation
1	Buttons	Control device functions (freely configurable).
2	Rocker switches	Control zoom and focus (horizontal and vertical configuration possible).
3	Joystick	Motorized XY fine adjustment (optional) in eight directions: 

3.6 Software description

3.6.1 General

All menus shown in these Instructions for Use are in English. The name and explanation of each function is specified in the respective language in the following lists.

Optional applications (licenses)

Various applications are available for the device which can be activated after obtaining a license. All licenses must be ordered separately.

Display and operation

Only activated applications are displayed on the touchscreen and can be configured and operated.

User-specific settings

User-specific settings are marked by a green dot placed in front of their designation.

These changes in settings are adopted only for the current user.

Please note the following: In order not to obstruct intraoperative changes of users, do not have the focus and zoom start values set user-specifically. Changes of the focus and zoom start values only become effective after the device is restarted.

3.6.2 Interactive buttons

All interactive buttons and their functions are described in the following table.

Symbol	Buttons	Explanation
	Live image	Touchscreen with live image, status bar and main menu bar
	Setup menu	All configurations and settings are accessible via this menu.
	switches	Function switched on
	switches	Function switched off
	Buttons	Triggers a function.
	Slide control	Continuously adjustable The values displayed (SI units) are rounded. They are intended to improve the display and do not represent a measurement function.

Symbol	Buttons	Explanation
	+/- buttons	Adjustable in steps
	Record button	Records and saves the current display.
	Video Recording start button	Video recording running
	Video Recording stop button	Video recording stopped
	Selection field	Selection between functions (blue: selected)
	functionalities	Selection of various functions (blue: selected)
	Loudspeaker on button	Switch loudspeaker on
	Loudspeaker off button	Switch loudspeaker off
	PIN lock deactivated	Unlock settings for handgrips and FCP.
	PIN lock activated	Lock Hand Grips and FCP Settings.
	Selection menu	Selection of various system languages
	functionalities	Create new user or patient
	functionalities	Edit patient/user profile
	functionalities	Filter patient data/filtering direction
	functionalities	Search for patient record
	functionalities	Import/Export data
	functionalities	Hide keyboard
	functionalities	Delete

Table 2: Interactive buttons

3.6.3 User interface of the KINEVO 900


Figure 29: Design of the KINEVO 900 user interface

1	Main menu (Home)	Navigation to the main menu with live image, status bar and main menu bar.
2	Edit Patient icon	Used to edit patient data.
3	Patient	Displays name of current patient, tapping on the button opens the patient directory.
4	edit a user	Used to edit user data.
5	Physician/user	Displays name of current user, tapping on the button opens the user directory.
6	Status Information	Displays current status information. The color of the hard disk or USB symbol changes from white to orange if only 20% of the storage capacity is available. Tapping on the button opens an overview [70].
7	Date, time	Displays current date and time.

Pos.	Name	Explanation
8	Setup menu	Opens the configuration menu.
9	Help	Opens the help window on the right-hand edge of the screen.
10	Light	Switches the illumination on / off.
11	Focus Light Link	Limits the maximum light intensity for the selected working distance.
12	Additional Light	Switch the additional light on/off.
13	Light Intensity	Slide switch for setting the light intensity.
14	Focus speed	Sets the speed of the motorized fine focusing.
15	Zoom speed	Set the zoom speed.
16	XY adjustment movement modes	Set the stand and speed movement modes.
17	Fluorescence setting	Select the activated fluorescence modules for BLUE 400, YELLOW 560 and INFRARED 800 with FLOW 800.
18	Position memory	Save position (bookmark) or move to stored position.
19	Number of photos	Displays the number of photos taken.
20	Image capture	Take still image (photo).
21	Number of videos	Displays the number of videos taken.
22	Duration of video recording	Shows the runtime of the current recording. The display is hidden if no video recording is active.
23	Video recording	Start and stop video recording.
24	Additional menu "Extended operation"	Offers additional adjustment options for the device.
25	Surgical microscope or endoscope	Tap on the [Surgical microscope] button to work only with the device. Tap on the [QEVO] button to switch to the QEVO Micro-Inspection Tool. [69]
26	AutoBalance	Tap on the [System] button to balance the device. Tap on the [Drape] button to perform fine balancing with a drape.
27	Video recording with low resolution	Start and stop video recording with low resolution. Suitable for recording long-term videos.

Pos.	Name	Explanation
28	AutoDrape	If the RFID label on your SMARTDRAPE has been scanned in, the system starts when the [AutoDrape] with drape suction button is tapped.
29	Foot control panel	Tap on the button in order to assign functions to the buttons of the foot control panel.
30	Start value of co-observer tube	Switching between face-to-face and lateral co-observer connection. Only available in the optional configuration for hybrid (optical-digital) visualization and in the case of active optical co-observation (see Tube [79]). In digital mode, a stereo glasses symbol refers to the digital setting.
31	Settings of current users	Displays the user-specific eyepiece settings. Only available in the optional configuration for hybrid (optical-digital) visualization and in the case of active optical co-observation (see Tube [79]). In digital mode, a stereo glasses symbol refers to the digital setting.
32	Hand Grips	Tap on the button in order to assign functions to the buttons of the hand grips.

All control elements, the status bar and the main menu bar can be hidden by tapping on the live image (full screen mode). Only the video runtime is displayed during a recording in full screen mode. If you tap on the monitor again, all control elements, the status bar and the main menu bar will reappear. The values displayed (SI units) are rounded. They serve to enhance the display and do not represent a measurement function.

3.6.4 QEVO user interface (option)


Figure 30: QEVO user interface

Poř.	Name	Explanation
1	Automatic Switch-off	Switch automatic deactivation on/off: ■ On (switch illuminated in blue): If the Micro-Inspection Tool is set aside and not moved, the camera and light source will switch off after approx. 40 seconds. To switch the camera and light source back on, press the [QEVO] button on the hand grip or the foot control panel. ■ Off (switch gray): The camera and light source of the Micro-Inspection Tool are permanently switched on. The presentation on the monitor and the MultiVision display (if activated) remain in QEVO mode.
2	White Balance	Perform white balance
3	Brightness	Set the image brightness on a scale from -5 to +5
4	QEVO Alignment	Adjust the camera image on the monitor to the actual viewing direction in the surgical field
5	Marker	Switch marker on/off: ■ On (switch blue): A marking in the image displays the direction of sight of the camera image. ■ Off (switch gray): No marking is displayed.

3.6.5 Status Information


Figure 31: Status information overview

Pos.	Symbol	Name	Explanation
1		Internal Hard Drive	Available internal hard disk space The color of the hard disk symbol changes from dark gray to orange if only 20% of the storage capacity is available.
2		USB	Available storage space on connected USB storage medium The color of the USB symbol changes from dark gray to orange if only 20% of the storage capacity is available.
3	Eject	Remove Hardware	The USB storage medium (e.g. USB stick or external Mini HDD USB 3.0 hard disk) may be removed. For reliable system functioning, this eject function must be used before each removal of a USB storage medium.
4		Connection status	Wireless FCP connected
5		Battery Status	Charging status of battery in FCP
6		Signal Strength	FCP signal strength
7		Hotspot	Hotspot connection (WiFi) ON/OFF

Pos.	Symbol	Name	Explanation
8		Show Hotspot Password	Hotspot password display ON/OFF
9		Web interface address	Displays the address of the web interface

3.6.6 "Patient data" menu


Figure 32: Design of "Patient data" menu

Pos.	Name	Explanation
1	Main menu (Home)	Navigation to the main menu with live image, status bar and main menu bar.
2	Patient	Displays name of current patient, tapping on the button opens the patient directory.
3	All	Displays all of the respective patient's recorded images, videos and data in the preview window.
4	Photos	Filters the preview window for photos; all of the respective patient's recorded images are displayed.
5	Videos	Filters the preview window for videos; all of the respective patient's recorded videos are displayed.
6	Data	Filters the preview window for data; all of the respective patient's recorded data are displayed (e.g. analysis data from the fluorescence mode INFRARED 800 with FLOW 800).

Pos.	Name	Explanation
7	Select All	Selects the entire contents of the respective preview window. This function is used e.g. to export or delete patient data.
8	Preview window display	Determines the size of the images in the preview window (large, medium, small).
9	Export	Exports all selected patient files.
10	Delete	Deletes all selected patient files.
11	Preview window	Displays all patient files as thumbnails.
12	Storage volume	Displays the total storage volume of all patient files in the footer bar.
13	Number of files	Displays the total number of all patient files in the footer bar.
14	Number of patients	Displays the total number of all patients in the footer bar.
15	Patient Directory	Displays the patient registry on the left side of the monitor.
16	Edit patient directory	Offers five options for working with the patient directory: Select all patients/Add new patients/Import patient files/Sort patients/Search for patients

3.6.7 Importing patient data from the DICOM Worklist

Main menu - Patient - Import  - DICOM


Figure 33: Input mask and DICOM Worklist display window

Pos.	Name	Explanation
1	Scheduled date	Period of time (selectable) for which the DICOM Worklist is to be called up.
2	Last name	Patient's last name.
3	First Name	Patient's first name.
4	Patient ID	Primary hospital identification number for the patient.
5	Accession Number	Department-specific case number which identifies the order for the "Imaging Service Request".
6	Req. Proc. ID	Requested Procedure ID: identification number that identifies the requested procedure in the "Imaging Service Request".
7	Modality	Device type with which the scheduled procedure step is to be performed.
8	AE title	AE title of the device with which the scheduled procedure step is to be performed.
9	Show Worklist	Tap on button to open the DICOM Worklist.

Pos.	Name	Explanation
10	Sch. Proc. Step Desc.	Scheduled Procedure Step Description: description specified by the hospital or classification of the scheduled procedure step.
11	Arrow button	Tap on button to display more details about the patient.
-	Req. Proc. Desc.	Requested Procedure Description: administrative description specified by the hospital or classification of the requested procedure.
-	Req. Proc. Code	Requested Procedure Code: code value that describes the requested procedure according to a specific coding scheme.
-	Ref. Phys. Name	Referring Physician's Name: name of the referring physician who requested the procedure.
-	Sch. Proc. Step Start	Scheduled Procedure Step Start: scheduled start time for the procedure step.
-	Sch. Protocol Code	Scheduled Protocol Code: code value that describes the scheduled protocol according to a specific coding scheme.

3.6.8 "Edit video recording" menu


Figure 34: "Edit video recording" menu

1	Video Recording	Displays the recorded video.
2	Fast forward	Used to run the video recording forward rapidly.
3	Cut out	Determines the end point of the timeline selection window, and therefore of the scheduled video recording.
4	Create New Video	Saves the cut video sequence as a new video.
5	Cut in	Determines the starting point of the timeline selection window, and therefore of the scheduled video recording/display.
6	"Video playback" slide control	The gray slide control keeps moving throughout the entire length of the video recording. Slide the gray slider to jump to an interesting area in the video recording without using the [Fast-forward] / [Rewind] buttons.
7	Save video sequence temporarily	Temporarily saves the currently cut video sequence and displays it as a preview.
8	Time marker	Marks the current position in the timeline selection window during video playback.
9	Preview window	Displays the cut video sequence(s).

Pos.	Name	Explanation
10	Timeline selection window	Displays the individual sequences/frames of the video. The start and end of the timeline selection window is bounded by blue bars. You can slide these blue bars to select and cut the sequence relevant for you in the video.
11	Start / Pause	Plays back the video recording. The playback of the video recording can be stopped by repeatedly tapping on the button.
12	Add Original Recording	Opens a preview window with all original video recordings of the patient.
13	Original Recording	Displays the original video recordings as a preview.

3.6.9 SmartLinks

In the "Settings Microscope" menu, four SmartLinks are available for your selection:

Figure 35: SmartLinks

Pos.	Name	Explanation
1	Focus Light Link	Limits the maximum light intensity for the selected working distance.
2	Focus Zoom Link	Automatically adapts the focusing speed to the magnification. The preselected focusing speed is automatically reduced when a higher magnification is used.
3	XY Zoom Link	Automatically adapts the XY travel speed to the magnification. The preselected travel speed is automatically reduced when a higher magnification is used.
4	Auto Brightness	Regulates the light intensity for constant brightness in the eyepiece depending on the working distance and magnification. Automatically adapts the luminous field diameter.

3.6.10 Light

In this menu, you configure the microscope illumination settings.


Figure 36: “Microscope Settings” menu, “Light” submenu

Pos.	Name	Explanation
1	Light	Switches the illumination on / off.
2	Additional Light	Switch the additional light on/off.
3	Lamp Service Life	Select which lamp is active (blue). Remaining operating time of the lamp until the next recommended lamp change (operating time per lamp = 500 hours).
4	Light intensity	Slide switch for setting the light intensity.
5	Light Warning Threshold	Set light warning threshold (factory setting: 25%). Display on monitor.
6	Change Light Intensity	Select speed of light intensity change (fast or slow).
7	SmartLinks: Focus Light Link	Limits the maximum light intensity for the selected working distance.
8	SmartLinks: Auto Brightness	Regulates the light intensity for constant brightness in the eyepiece depending on the working distance and magnification. Automatically adapts the luminous field diameter.

3.6.11 Tube

This menu is only available in the optional configuration of the microscope for the hybrid (optical-digital) visualization.

Activate the optical co-observation and select the tubes and eyepieces used for the surgeon. The total magnification is calculated with these stored values and displayed on the monitor.

Figure 37: "Microscope Settings" menu, "Tube" submenu

Pos.	Name	Explanation
1	System Mode	Visualization mode selection: ■ Digital: purely digital visualization - optical co-observation not possible ■ Optical: hybrid visualization - co-observation through tubes and on the monitor Only available in the optional configuration for hybrid (optical-digital) visualization. The selection also influences options and displays in other menus.
2	Focal Length	Focal length of tube in use: (170 mm for tiltable tube and foldable tube, 260 mm for foldable tube with magnification).
3	Eyepiece magnification	Magnification of the eyepieces used: 10x or 12.5x.
4	Magnification Changer	Magnification of the optional 3-position magnification changer.
5	Co-observer (start value)	Selection of co-observer: ■ Lateral image outputs: Right and left for stereo co-observation module or photo adapter for external camera ■ Opposite image outputs: Face-to-face for tiltable tube or foldable tube

Pos.	Name	Explanation
		The “co-observer (start value)” can be adjusted on the monitor in the additional menu “Extended operation” and manually on the microscope body.

3.6.12 Focus

In this menu, you configure the focus settings.


Figure 38: “Microscope Settings” menu, “Focus” submenu

Pos.	Name	Explanation
1	Start value	The microscope starts with this working distance after it is switched on.
2	Speed	Set the speed of the motorized fine focus adjustment.
3	Autofocus (AF) Speed	(option) Adjust the autofocus speed.
4	AF Auto Activation (brakes)	(Option) On: Autofocus is automatically activated each time the brakes are closed. Off: Autofocus only via configurable button on the handgrip or foot control panel. Autofocus is not available during a few specific applications. A corresponding acoustic signal is then emitted.
5	Laser Focus Aid	(option) Laser Focus Aid on: The focusing laser is switched on when the brakes are released or the motorized fine focus adjustment is actuated. Laser Focus Aid off: The focusing laser is switched off.
6	Focus Lock	For the application of an optional micromanipulator with a fixed working distance.

Pos.	Name	Explanation
		Focus Lock on: Autofocus inoperative. Focus rocker switch on handgrip/foot switch out of operation. Focus Lock off: Autofocus functioning. Focus rocker switch on handgrip/foot control panel active.
7	SmartLinks: Focus Zoom Link	Automatic adaptation of the focusing speed to the magnification. The preselected focusing speed is automatically reduced when a higher magnification is used.

3.6.13 Zoom

In this menu, you configure the zoom settings.

Figure 39: "Microscope Settings" menu, "Zoom" submenu

Pos.	Name	Explanation
1	Gamma (start value)	Set the zoom factor (0.4 ... 2.4) at which the zoom system is to start after the system is switched on.
2	Speed	Set the zoom speed.
3	SmartLinks: XY Zoom Link	Automatic adaptation of the motorized XY travel speed to the magnification. The preselected travel speed is automatically reduced when a higher magnification is used.

3.6.14 Hand Grips

In this menu, you configure the controls of the handgrips.

Figure 40: "Operation settings" menu, "Handgrips" submenu

Pos.	Name	Explanation
1	Button [default factory setting]	Tapping this button resets the programming of all handgrip functions to the factory settings.
2	Selection window	Appears after a handgrip button is tapped: Displays all available functions for each handgrip button. Tap a function to apply it for the selected button.
3	Enter PIN	By activating this button, you can protect your button assignment for the handgrips with a four-digit PIN.
A, B, C, D	Configurable buttons	Tapping on the respective button opens a selection window on the right side where the functions available to you can be selected.
S	XY Joystick	For motorized XY fine adjustment: The XY joystick of the right handgrip can have certain special functions, depending on the option activated. Messages displayed on the touchscreen can be acknowledged by pressing the center position of the joystick button on the left handgrip.
E, F, G, H	Rocker switch for focus and zoom	Tapping on the buttons opens a selection window on the right side for selecting the direction (+/-) or changing the zoom/focus functions.

Pos.	Name	Explanation
SB	[SB] ("Selected Brakes") button	Tapping on the button opens, on the right-hand side, a selection window for programing the button with the brake function for the stand brakes, the brake function for the microscope brakes or the PointLock function.
AB	[AB] ("All Brakes") button	The button is not configurable. As long as the button is pressed, the entire system can move freely. When you let go of the button, the system locks stably and with vibration damping.

3.6.15 Foot Control Panel (FCP)

In this menu, you configure the controls of the FCP (foot control panel).


Figure 41: "Operation settings" menu, "FCP" submenu

Pos.	Name	Explanation
1	Button [default factory setting]	Tapping this button resets the programming of all handgrip buttons to the factory settings.
2	Selection window	Appears after a foot control panel button is tapped: Displays all available functions for each foot control panel button. Tap a function to apply it for the selected button.
3	Enter PIN	By activating this button, you can protect your button assignment for the foot control panel with a four-digit PIN.
S	XY Joystick	Motorized XY fine adjustment.
A, B, C, D, E, F	Configurable buttons	Tapping on the respective button opens a selection window on the right side where the functions available to you can be selected.
G, H, I, J	Rocker switch for focus and zoom	Tapping on the buttons opens a selection window on the right side for selecting the direction (+/-) or changing the zoom/focus functions.

3.6.16 Rock Switch

In this menu, you configure the controls of the rocker foot switch.


Figure 42: "Operation settings" menu, "Rocker foot switch" submenu

Pos.	Name	Explanation
1	Button [default factory setting]	Tapping this button resets the programming of the rocker foot switch to the factory settings.
2	Selection window	Appears after rocker foot switch is tapped: Displays all available functions for each rocker foot switch button. Tap a function to apply it for the selected button. When the following functions are selected for a rocker foot switch, the second rocker foot switch is automatically assigned: ■ Light: Light brighter (top button) / darker (bottom button) ■ Zoom: + (top button) / - (bottom button) ■ Focus: + (top button) / - (bottom button)
3	Enter PIN	By activating this button, you can protect your button assignment for the rocker foot switch with a four-digit PIN.

3.6.17 Suspension system

In this menu, you configure the movement modes, the park position and the drape position.


Figure 43: "Settings Stand" menu

Pos.	Name	Explanation
1	XY Mode	Selection of the movement modes for motorized fine adjustment using the XY joystick on the handgrip or on the foot control panel: ■ Microscope: Microscope movement in three axes ■ Stand: XY stand movement in the focal point plane ■ PointLock: adjustment around the focal point The setting can also be carried out in the main menu.
2	Start Park Position	Device moves to park position. To do this, the right button of the XY joystick on the hand grip or the foot control panel must also be kept pressed.
3	Start Drape Position	Device moves to drape position. To do this, the right button of the XY joystick on the hand grip or the foot control panel must also be kept pressed.
4	Save Current Position as Drape Position	Save drape position.

Pos.	Name	Explanation
5	Microscope	Preselect speed for movements in this XY mode.
6	Suspension system	Preselect speed for movements in this XY mode.
7	PointLock	Preselect speed for movements in this XY mode.
8	PositionMemory	Preselect speed for movements in the PositionMemory function.

3.6.18 Displays

In this menu, you configure the settings of the video outputs.


Figure 44: "Display settings" menu

Pos.	Name	Explanation
1	Video format for second system monitor and optional external monitor	Video format (only available with integrated 3D video system / stereo option) ■ 2D ■ 3D Output ■ Copy System Monitor ■ Application-specific
2	Output HD-SDI	■ Progressive (full image) ■ Interlaced
3	DisplayPort output	■ Left-hand camera ■ Right-hand camera
4	Selection window	Appears after tapping on video output: Shows all available output options for each output. Tap an option to apply it for the selected output.

3.6.19 Audio

You can configure the audio settings in this menu.


Figure 45: “Settings Display” menu

Pos.	Name	Explanation
1	Internal Micro- phone: Use for Video Recording	Switch microphone on/off.
2	Internal Micro- phone: Forward to Line-out	Switch microphone on/off.
3	External Input (line- in): Use for Video Recording	Switch microphone on/off. To avoid disruptive noises, the internal microphone should be switched off while an external microphone is being used.
4	External Input (line- in): Forward to Line-out	Switch microphone on/off.
5	External Input (line- in): Forward to Internal Speaker	Switch microphone on/off. Adjust volume.
6	Volume: Video Playback	Switch playback on/off. Adjust volume.
7	Volume: System Sounds	Switch playback on/off. Adjust volume.

3.6.20 Photo

In this menu, you configure the photo settings.


Figure 46: "Settings Photo" menu

Pos.	Name	Explanation
1	Record	Selection: - Single images with internal video camera - Images with external photo camera (DSLR, option)
2	Storage Format (internal camera)	Storage format for images taken: - PNG format - JPEG format
3	Light Intensity for External Camera	Utilization of maximum xenon light output (100%) for taking photos.

3.6.21 Video Camera

In this "White Light" menu, you configure the settings of the video camera. You can optionally configure the video settings for BLUE 400 and YELLOW 560.

The menu setting options vary depending on whether the camera shutter control is set at "Auto" or "Manual".

Options in automatic shutter control


Figure 47: "Video Camera settings" menu, automatic shutter control

Pos.	Name	Explanation
1	Preview window	Video image for checking the settings.
2	Shutter Control	■ Auto: Automatic control of the exposure time ■ Manual: Manual control of the exposure time
3	Exposure Control	■ Spot ■ Full
4	Auto Brightness	Set the brightness level of the automatic brightness control.
5	Color Saturation	Set the color saturation.
6	Hue	Set the hue.
7	Black Level Offset	With a higher black level, dark image details are displayed brighter. This can lead to an enhanced recognizability.
8	Red Fine Adjustment	Set the red value.
9	Blue Fine Adjustment	Set the blue value.
10	White Balance	Start white balance.
11	Shutter Control BLUE 400	■ Auto: Automatic control of the exposure time ■ Manual: Manual control of the exposure time
12	Auto Brightness BLUE 400	Set the brightness level of the automatic brightness control for BLUE 400.
13	Shutter Control YELLOW 560	■ Auto: Automatic control of the exposure time ■ Manual: Manual control of the exposure time
14	Auto Brightness YELLOW 560	Set the brightness level of the automatic brightness control for YELLOW 560.

Options in manual shutter control


Figure 48: "Video Camera" menu, manual shutter control

Pos.	Name	Explanation
1	Preview window	Video image for checking the settings.
2	Shutter Control	■ Auto: Automatic control of the exposure time ■ Manual: Manual control of the exposure time
3	Manual Shutter	Set exposure time manually
4	Color Saturation	Set the color saturation.
5	Hue	Set the hue.
6	Black Level Offset	With a higher black level, dark image details are displayed brighter. This can lead to an enhanced recognizability.

Pos.	Name	Explanation
7	Red Fine Adjustment	Set the red value.
8	Blue Fine Adjustment	Set the blue value.
9	White Balance	Start white balance.
10	Shutter Control BLUE 400	■ Auto: Automatic control of the exposure time ■ Manual: Manual control of the exposure time
11	Manual Shutter BLUE 400	Set exposure time for BLUE 400 manually.
12	Shutter Control YELLOW 560	■ Auto: Automatic control of the exposure time ■ Manual: Manual control of the exposure time
13	Manual Shutter YELLOW 560	Set the desired exposure time for YELLOW 560 manually.

3.6.22 Connections

You can inspect or configure the connection settings in this menu. The settings in this menu are password-protected and can only be changed by authorized users with IT-Admin system rights [149]. For other users the settings cannot be edited.

"Computer Name" submenu


Figure 49: "Connectivity Settings" menu, "Computer Name" submenu

Item	Name	Description
1	Change computer name (restart required)	Change computer name, changes will take effect after restart
2	[Save] button	Apply setting.

"Network Drive" submenu

Figure 50: "Connectivity Settings" menu, "Network Drive" submenu

Item	Name	Description
1	Connect to network drive	Activate connection to the network drive.
2	Host	Enter host name.
3	User Name	Enter user name.
4	Password	Enter the password.
5	Workgroup / Domain	Enter workgroup.
6	Share Name	Enter share name (sys. admin data).
7	[Test Connection] button	Test connection to network drive.
8	[Save] button	Save entry.

"Ping IP Address" submenu


Figure 51: "Connectivity Settings" menu, "Ping IP Address" submenu

Item	Name	Description
1	IP Address	Enter IP address for the Ping test.
2	[Ping] button	Start ping test.

"LAN" submenu

Figure 52: "Connectivity Settings" menu, "LAN" submenu

Item	Name	Description
1	Connection status	Status of the LAN connection
2	Connect to LAN	Activate LAN connection.
3	MAC Address	Display the MAC address
4	DHCP	Configure network connection via DHCP.
5	IP Address	Enter IP address for the LAN connection. NOTE! Do not enter an IP address from the following areas: ■ 192.168.173.xxx (is used for the ZEISS KINEVO 900 hotspot) ■ 192.0.2.xxx (is used for an internal network connection with QEVO (Option))
6	Subnet Mask	Enter subnet mask for the LAN connection.
7	Standard Gateway	Enter Standard Gateway for the LAN connection.
8	Use DNS Server Address Automatically	The DNS server address is automatically activated with DHCP.
9	Primary DNS server	Address of the primary DNS server
10	Secondary DNS Server	Address of the secondary DNS server
11	[Save] button	Apply settings.

"WLAN" submenu (WPA2/PSK)

The "WLAN" submenu varies depending on which authentication log is selected.



Figure 53: "Connectivity Settings" menu, "WLAN" submenu with WPA2/PSK authentication

Item	Name	Description
1	Connection status	Status of the WLAN connection
2	Connect to WLAN	Activate WLAN connection.
3	MAC Address	Display the MAC address
4	Authentication	Select authentication standard (here: WPA2/PSK).
5	SSID	SSID / enter network name.
6	Password	Enter network password.
7	DHCP	Configure network connection via DHCP.
8	IP Address	Enter IP address for the WLAN connection. NOTE! Do not enter an IP address from the following areas: ■ 192.168.173.xxx (is used for the ZEISS KINEVO 900 hotspot) ■ 192.0.2.xxx (is used for an internal network connection with QEVO (Option))
9	Subnet Mask	Enter subnet mask for the WLAN connection.
10	Standard Gateway	Specify Standard Gateway for the WLAN connection.

Item	Name	Description
11	Use DNS Server Address Automatically	The DNS server address is automatically activated with DHCP.
12	Primary DNS server	Address of the primary DNS server
13	Secondary DNS Server	Address of the secondary DNS server
14	[Save] button	Apply settings.

"WLAN" submenu (WPA2/PEAP)

The "WLAN" submenu varies depending on which authentication log is selected.

Item	Name	Description
1	Connection status	Status of the WLAN connection
2	Connect to WLAN	Activate WLAN connection.
3	MAC Address	Display the MAC address
4	Authentication	Select authentication standard (here: WPA2/PEAP).
5	SSID	SSID / enter network name.
6	User Name	Enter user name.
7	Password	Enter password.
8	[Import certificate] button	Load WLAN certificate from USB storage medium.
9	DHCP	Assign a network connection via DHCP.

Item	Name	Description
10	IP Address	Enter IP address for the WLAN connection. NOTE! Do not enter an IP address from the following areas: ■ 192.168.173.xxx (is used for the ZEISS KINEVO 900 hotspot) ■ 192.0.2.xxx (is used for an internal network connection withQEVO(Optional))
11	Subnet Mask	Enter subnet mask for the WLAN connection.
12	Standard Gateway	Specify Standard Gateway for the WLAN connection.
13	Use DNS Server Address Automatically	The DNS server address is automatically activated with DHCP.
14	Primary DNS server	Address of the primary DNS server
15	Secondary DNS Server	Address of the secondary DNS server
16	[Save] button	Apply settings.

"Hotspot" submenu

Figure 54: "Connectivity Settings" menu, "Hotspot" submenu

Item	Name	Description
1	Enable Hotspot	Switch hotspot on / off. WLAN must be activated to switch on the hotspot.
2	SSID	Define SSID / network name.
3	Password	Define network password.
4	[Save] button	Apply settings.

"Web Interface" submenu


Figure 55: "Connectivity Settings" menu, "Web Interface" submenu

Item	Name	Description
1	Password	Define Web interface password. Factory setting: "ZEISS".
2	[Save] button	Apply settings. For the connection to the web interface, WLAN or LAN must be activated.

"DICOM" submenu


Figure 56: "Connectivity Settings" menu, "DICOM" submenu

Item	Name	Description
1	Station Name	Enter recognizable device names.
2	AE Title	Enter DICOM names of the AE (Application Entity).
3	Port	Port which is used for the connection with DICOM (pre-set)
4	Non-Listening Mode	Block incoming DICOM connections. Note: When option activated, PACS does not respond to DICOM C-ECHO.
5	MWL AE Title	Enter work list names.
6	MWL Hostname	Enter host name.
7	MWL Port	(pre-set)
8	Storage AE title	Enter AE title of DICOM storage.
9	Storage Hostname	Enter DICOM Storage hostname.
10	Storage port	(pre-set)
11	Institution Name	Enter name of clinic/institution.

Item	Name	Description
12	[Test Connection] button	Test DICOM connection.
13	[Save] button	Apply settings.

"Advanced DICOM Settings" submenu

Item	Name	Description
1	Association Idle Timeout [s]	maximum period of time in which the DICOM connection may remain unused without being separated
2	DIMSE Response Timeout [s]	maximum period of time in which (with active DICOM connection) no DIMSE response from DICOM must be received
3	Connect Timeout [s]	maximum period of time until a DICOM connection is established
4	Log Level	Setting the log level for the used components of the DICOM network
5	Modality	List of modality codes / modality types that may be selected as modality [223]
6	Max Query Results	maximum number of results to a query
7	Transfer Syntax 2D	MPEG transfer syntax for 2D video files
8	Transfer Syntax 3D	MPEG transfer syntax for 3D video files
9	Char Set	Character set which is used for the DICOM connection
10	[Save] button	Apply settings.

"ZEISS Smart Services" submenu

Figure 57: "Connectivity Settings" menu, "ZEISS Smart Services" submenu

Item	Name	Description
1	Connection status	Status of the ZEISS Smart Services connection
2	Activate	Activate ZEISS Smart Services.
3	Automatic SCU Detection	Find SCU server automatically.
4	IP Address	Enter IP address of the SCU server.
5	Port	Enter port for communication with the SCU server.
6	[Save] button	Apply settings.

3.6.23 Video recording

In this menu, you configure the settings for video recordings.


Figure 58: "Settings Video Recording" menu

Pos.	Name	Explanation
1	Parallel Recording	The data are always stored on the internal hard drive. In addition, parallel storage on a connected USB storage medium or a configured network drive (USB or network drive) also is possible.
2	3D Recording	Option for a stereo video camera. Selection whether the internal recording is to be made in 3D format (ON).
3	SmartRecording	If "SmartRecording" is activated, video data of the operation are recorded in a buffer which is automatically cleared and restarted every five minutes. You can start and stop the recording of the video data from the buffer at any time.
4	Video Streaming	Video recordings can be live streamed to a preset web address or device via LAN or WLAN.
5	Video Streaming	Selection of video resolution. ■ HD Video: High resolution (1920 x 1080) ■ LowRes Video: Low resolution (960 x 540)
6	Web Address for Video Streaming	Entry of intended web address for video streaming.

3.6.24 Fluorescence

(Option) Fluorescence settings for BLUE 400, YELLOW 560 and INFRARED 800 with FLOW 800.

Only activated fluorescence options are displayed on the touch-screen and can be configured.


Figure 59: "Settings Fluorescence" menu

Pos.	Name	Explanation
1	Fluorescence mode	Selection of the fluorescence mode which is activated with the configured hand grip button ■ BLUE 400 ■ YELLOW 560

Pos.	Name	Explanation
		■ IR 800 ■ IR 800 with FLOW 800 By pressing on the touchscreen, the configured handgrip button or the FCP button, you can activate the respective fluorescence application.
2	Light Intensity Start Value BLUE 400	Selecting the light intensity for BLUE 400
3	Light Intensity Start Value YELLOW 560	Selecting the light intensity for YELLOW 560
4	AutoGain	Manual or automatic control of camera gain ■ Auto: The camera gain is controlled automatically. ■ Manual: The camera gain is controlled manually.
5	AutoDetection	AutoDetection detects the influx of the fluorescence dye during the recording. ■ ON: The black leader without any fluorescence signal is hidden during the replay. ■ OFF: The replay is played back incl. the black leader.
6	Skip Setup	When the option is activated, INFRARED 800 starts immediately in the recording mode, the setup phase is skipped.
7	AutoZoom	On: The total magnification is automatically changed to the preconfigured value without a notice dialog. Off: A setup dialog for changing the total magnification outside of the preconfigured range is displayed.
8	AutoZoom (start value)	Magnification value to which the microscope is set for the start of the recording of the INFRARED 800 video.
9	Gain (start value)	Set the gain start value for a manual camera gain.
10	Short Replay	The short replay is played back repeatedly until the preset number of short replays has been reached or the button that was configured for the fluorescence application is pressed on the handgrip or the foot control panel. The number of short replays is automatically limited to max. 25.

Pos.	Name	Explanation
11	Short Replay (duration)	The length of the short replay is adjustable from 6 to 30 seconds. The adjusted length defines the video sequence used for the FLOW 800 card.
12	Long Replay	The recorded fluorescence video is played back in full length until the preset number of replays has been reached or the button that was configured for the fluorescence application is pressed on the handgrip or the foot control panel.
13	MultiVision Overlay	Live display of the infrared image (during recording) ■ On: Display on the monitor and as MultiVision Data Injection ■ Off: Display on the monitor only
14	Display on Monitor	Live display of the infrared image on the touchscreen (during recording) ■ Full Image ■ Overlay ■ Picture in Picture
15	MultiVision Data Injection	Show FLOW 800 evaluation

3.6.25 MultiVision

In this menu, you configure the settings for MultiVision (integrated data injection, option).


Figure 60: "Settings MultiVision" menu (option)

Pos.	Name	Explanation
1	MultiVision	Switch the MultiVision function (option) on/off
2	Source	Selection of MultiVision data source for display: ■ None ■ Navigation: Data from the navigation system are displayed. ■ QEVO / Video: Video recording of the Micro-Inspection Tool or the integrated video camera is shown.
3	Status Displays	Selection of displays which can be shown in the MultiVision display with a corresponding symbol: ■ Focus: current working distance ■ Zoom: current magnification ■ Light intensity: current light intensity ■ Reticle: marking of the current focus point

Pos.	Name	Explanation
		■ Recording: REC symbol ■ Microphone: current status of microphone ■ XY mode: current movement mode of XY adjustment ■ Camera (16:9): Marking of the camera recording area
4	Brightness	Set the brightness of the data injection
5	Show MultiVision Demo Videos	Selection of different MultiVision demo videos. These MultiVision demo videos can be displayed with the shutter open or closed.
6	QEVO Automatic Switch-off (Auto rest mode)	Switch automatic switch-off for QEVO on/off
7	QEVO Marker	Switch QEVO marker on/off
8	Plug Cycles	Shows how often a QEVO has already been plugged in.
9	Video Source	Select video source (output 5): ■ QEVO ■ External
10	Navigation (standard)	The navigation system is connected directly to the device: Activate / deactivate connection to navigation system with standard functions
11	Navigation (extended)	The navigation system is connected to the device via the network: Activate / deactivate connection to navigation system with extended functions. As soon as the connection to the navigation system has been activated, the current IP address appears in the "Navigation (extended)" field.

3.6.26 Info

This menu displays a variety of information regarding the system (serial numbers, software versions, license information).



Figure 61: "Settings Info" menu

Pos.	Name	Explanation
1	Serial Number	Displays serial number of device.
2	Name / Version	Displays the following values: ■ Product name ■ Software version
3	Permanent Licenses	Lists all licenses installed in the device.
4	Export License Information to USB	Tap on the [Start] button to export licenses to a USB storage medium.
5	Radio Frequency Approval	Tap on the [Open] button to display all radio frequency approvals available on the device. Note: You can also request the latest version of the radio frequency approval from your ZEISS representative.
6	Windows End User License Agreement (EULA)	Tap on the [Open] button to display the Windows End-user license agreement (English).

3.6.27 Extras

In this menu, you can perform license imports and log file exports.


Figure 62: "Settings Extras" menu

Pos.	Name	Explanation
1	License import	Import additional licenses
2	Export Log Files	When you press this button, the log files are exported to a connected USB storage device. You can then send this data to ZEISS Service as an email attachment.
3	Backup on Internal Hard Drive	Data is saved on internal hard drive
4	Most Recent Backup from:	Date of the last backup
5	Diagnostic Package	Send diagnostic data to ZEISS Note: This function is only available if ZEISS Smart Services are activated [241].
6	Remote Service Access	Enables a direct data exchange between the device and ZEISS Service Note: This function is only available if ZEISS Smart Services are activated [241].

3.6.28 Service PC

 
The settings in this menu are password-protected and may only be seen and changed by authorized users who have an IT Admin password [149].

3.6.28.1 Service PC / Overview

 Service PC - Overview

Figure 63: "Settings Service PC" menu, "Overview" tab

Pos.	Name	Explanation
1	Software information	Information on system software
2	Firmware information	Information on system component firmware

3.6.28.2 Service PC / Service

Service PC - Service


Figure 64: "Settings Service PC" menu, "Service" tab

Pos.	Name	Explanation
1	Synchronize Patient Data	Synchronize patient data.
2	Export User	Export the existing user profile.
3	Import User	Import the existing user profile.
4	Maintenance (date)	Date of last maintenance call.

3.6.28.3 Service PC / Settings

Service PC - Settings


Figure 65: "Settings Service PC" menu, "Settings" tab

Pos.	Name	Explanation
1	Activate User Password	Activate / deactivate user password (default setting: activated).
2	Standard Password	Change standard password for all users.

Pos.	Name	Explanation
3	System Language	Select System Language. The system language is applied before the successful user login and in the password-protected setting area.
4	Video Frequency	Set to 50 Hz or 60 Hz.
5	Date and Time	Set and [Accept].
6	Time Format	Select: 24h or 12h.
7	Date Format	Select: DD.MM.YYYY or MM/DD/YYYY.

3.6.28.4 Service PC / Log files

⚙️ Service PC - Log files


Figure 66: "Settings Service PC" menu, "Log Files" tab

Pos.	Name	Explanation
1	Export Log Files	Export log files to a USB storage device.
2	Export Audit Files	Export audit files to a USB storage medium.
3	[Load Log File]	Load log files from USB storage medium.
4	File type	Select type of log files to be loaded.

3.6.28.5 Service PC / Configuration

⚙ Service PC - Configuration


Figure 67: "Settings Service PC" menu, "Configuration" tab

Pos.	Name	Explanation
1	Delete All Licenses	Delete all licenses installed on the device.
2	Reimport Licenses	Reimport previously exported licenses.

Empty page, for your notes

4 Installation

4.1 Requirements

Installation requires special knowledge and skills.

An installation by persons not authorized by ZEISS can lead to the injury of patients and operators, as well as to property damage.

Action

- Have the installation and initial startup performed only by persons authorized by ZEISS.

CAUTION!

Risk of injury caused by uncontrolled device movements.

The device may make uncontrolled swinging movements if the suspension arm and the surgical microscope are not in the transport position.

- Place the device in the transport position before moving it.

4.2 Attaching the tubes and eyepieces

The information in this paragraph only applies to the optional configuration of the microscope as a system for hybrid (optical-digital) visualization.


Figure 68: Attaching the tubes and eyepieces

1	Wide-field eyepieces	2	Tiltable tube, optional foldable tube
3	Angle optics (spine adapter)*	4	Fastening screw, magnification changer
5	Magnification changer, 3-position*	6	Dust protection cover
7	Straight tube	8	Fastening screw
9	Transport protection cover	10	Fastening screw, angle optics

* When assembling the microscope, be careful not to mount the 3-position magnification changer (item 5) and the angle optics (item 3) on the microscope body at the same time. In this case, the microscope body would be too heavy for the autobalance function.

NOTE

Wide-angle eyepieces with magnetic coupling!

Please keep in mind the usual rules for the handling of magnets with regard to eyepieces that have been removed from the tube:

- ▶ Do not place the eyepiece near instruments which may be magnetizable.
- ▶ Do not place the eyepiece on sensible devices such as infusion pumps, pacemakers, measuring devices or magnetic data media such as floppy disks, audio and video tapes or credit cards.
- ▶ Always store the eyepiece in its original packaging when not using it.

⚠ CAUTION!

Never stare directly into the sun or a light source with the tube!

An excessive light intensity can damage the retina of the eye.

- ▶ Never look at the sun or a light source through the eyepieces or the binocular tube.

⚠ CAUTION!

Risk of injury to the patient caused by lowering of the surgical microscope or falling parts!

Never change modules or permissible components during a surgical procedure or above the patient!

- ▶ Always make sure that the maximum permissible load is not exceeded
- ▶ Balance the completely equipped surgical microscope before each application.

Action

1. Unscrew the fastening screw on the microscope body several turns.
2. Remove the dust protection cover.
3. If necessary: Place the angle optics with dovetail (spine adapter) or the magnification changer (option) on the microscope body and tighten the fastening screw on the microscope body hand-tight.
4. Mount the binocular tube and tighten the fastening screw on the microscope body or on the angle optics with dovetail (spine adapter) hand-tight.
5. Insert widefield eyepieces as far as they will go in the mounts provided on the binocular tube.
6. Remove the transport protection cover from the objective lens.
7. Enter the tube focal lengths and eyepiece magnifications used in the "Settings Surgical Microscope" menu so that the device can calculate the correct total magnification.

4.3 Using an integrated 3D video system for observation without eyepieces (option)

The fully digital configuration of ZEISS KINEVO 900 enables exclusively digital operation with an external 3D-video monitor. The optical co-observation through tubes and eyepieces cannot be used in this configuration. In the optional configuration as system for hybrid (optical digital) visualization, you have the possibility of using your ZEISS KINEVO 900 in the digital system mode for a complete observation without an eyepiece.

In this case, the surgical field is displayed on the second 3D video monitor or on a separate, external 3D video monitor via the 3D video system. With this type of application, the magnification, and therefore the resolution on the monitor, is lower than it is when using the surgical microscope with tube and eyepieces.

CAUTION!

Do not use images and videos for diagnostic purposes!

The monitors are neither calibrated nor designed for diagnostic purposes.

- The video contents and images displayed on the monitors (cut sequences) must not be used for diagnostic purposes. They may contain deviations with respect to scale, shape, contrast and color.

To use the surgical microscope - which is configured for optional hybrid visualization - as a digital 3D video system for observation without eyepieces, proceed as follows:

Prerequisite

- ☒ The microscope is configured as a system for hybrid (optical-digital) visualization.

Action

1. Unscrew the fastening screw on the microscope body several turns.
2. Remove the binocular tube with the widefield eyepieces from the microscope body or from the angle optics with dovetail (spine adapter).
3. Remove the angle optics with dovetail (spine adapter) or the magnification changer (option) from the microscope body.
4. Place the Digiskop cover on the microscope body and tighten the fastening screw hand-tight.
5. Unscrew the fastening screw of the binocular assistant's tube several turns.
6. Remove the binocular assistant's tube with the widefield eyepieces from the microscope body.
7. Place the assistant's Digiskop cover on the microscope body and tighten the fastening screw hand-tight.
8. Switch the device on and log into the system software.
9. Tap on  Settings →  Microscope → Tube [79].

10. Select the "Digital" option in the "System mode" field.
11. Leave "Settings" and check the image and color quality of the video display on the 3D video monitor.
 - ⇒ If the image and color quality are perfect, you can start your operation.

 WARNING!

Failure of digital image

If the 3D video system fails and the image is no longer displayed on the 3D monitor, the 3D video system may no longer be used for the operation.

- ▶ Complete the operation with a backup device or with an alternative means of visualization (fully digital configuration: loupe; hybrid configuration: tube).
- ▶ Unplug the device from the electric power supply and contact ZEISS Service.

4.4 Attaching the documentation/co-observer facility

The information in this paragraph only applies to the optional configuration of the microscope as a system for hybrid (optical-digital) visualization.


Figure 69: Attaching the documentation/co-observer facility

1	[Pivoting mirror] adjustment knob	2	Co-observation assembly e.g. photo adapter for DSLR or stereo co-observer
3	Dust cover	4	Knurled ring, right Opposite knurled ring on the left not shown

Action

1. Loosen the knurled ring concerned.
2. Remove the dust protection cap.
3. Push the co-observation module into the mount of the image output as far as it will go.
The mount of the image exit port is fitted with guide projections.
4. Screw the knurled ring onto the co-observation module and tighten the knurled ring securely.
5. Set the [pivoting mirror] adjustment knob on the microscope body for the lateral co-observer connection ②.
6. Or configure the pivoting mirror setting on the monitor:
Tap on ⚙ Settings → ⏏ Microscope → Tube [79].
7. In the [Co-observer (start value)] field, activate the function "Lateral".

4.5 Mounting the mouth switch

The information in this paragraph only applies to the optional configuration of the microscope as a system for hybrid (optical-digital) visualization.


Figure 70: Mounting the mouth switch

1	Protective cap, 3 pcs	2	Fastening screw, 3 pcs
3	Mouth switch for 180° tiltable tube and foldable tube	4	Mouth switch cable with connector
5	Angle adapter for mouth switch		

Action

1. Remove the three protective caps from the tube.
⇒ Threads for attaching the mouth switch are located underneath the protective caps.
2. Position the mouth switch on the tube and attach it with three screws (included in the scope of supply).
3. Tighten the three screws hand-tight using a 4 mm Allen key.
4. Insert the mouth switch cable in the angle adapter for the mouth switch.
5. Connect the angle adapter for the mouth switch to the connection socket on the microscope.
6. Set the height, tilt and distance of the mouth switch via the relevant locking screws and properly tighten these screws.
7. Perform an autobalance for the system.
8. Press one of the brake release buttons on the back of the hand grips to activate the mouth switch.

Note:

Each time the device is rebooted, the mouth switch will be activated only after one of the brake release buttons on the back of the hand grips is pressed.

You will find an exact description of the mouth switch in the Instructions for Use of the Mouth Switch G-30-1469.

4.6 Moving the device


Figure 71: Moving the device

1	Transport direction 	2	Locking mechanism tab 
3	Pedal for straight-ahead movement 	4	Holder for cable and foot control panel
5	Transport handles	6	Parking position [190]

CAUTION!

Risk of crushing!

Fingers may be crushed between the vertical arm and the horizontal arm.

- ▶ Never touch the area between the vertical arm and the horizontal arm while moving the device.
- ▶ Use the transport handles to move, push and position the device.

CAUTION!

Risk of toppling!

The device may topple over and injure persons when being moved over thresholds and obstacles.

- ▶ Be careful when pushing or pulling the device over thresholds.
- ▶ Push or pull the device **grasping it only by its transport handles** in the transport direction.
- ▶ Be extremely careful when moving over inclines.
- ▶ Do not park the device on sloped surfaces.

Action

1. Tap on  Settings →  Stand → and, in the -Move to park position field, → the -Start button
2. To activate the travel, press the XY joystick button on the right handgrip to the right  and hold it down until the park position has been reached. The joystick on the FCP also can be used to activate the travel.
3. The device now automatically travels to the park position.

⇒ When the park position has been reached, this is acknowledged by a signal tone.

4. Remove the power plug from the wall socket.
5. Hang the foot control panel on the holder provided for this purpose.
6. Hang the power cord on the holder provided for this purpose.
7. Release the locking button and push the device in the direction of the arrow (= transport direction).
8. Press the straight-ahead travel button:
When steering the device over longer distances and in a straight-ahead direction.
9. When the locking button is pressed gently, all four castors again become steerable, e.g. in order to position the device at its usage site.
10. Press the locking button:
For locking at the usage site or in a park position.

4.7 Connecting the device to the power supply

A yellow locking mechanism is attached to the power cord to prevent it from being accidentally disconnected.

Power cords in The Netherlands and Switzerland:

The power cords used in The Netherlands and Switzerland (order numbers: 000000-0603-410 and 000000-0584-947), do not have a locking mechanism. A potential equalization cable with a connector is attached to both ends of the power cord in each case.

CAUTION!

Danger! Electrical voltage!

When the power cord is connected to the power outlet, the device is in the Standby mode. The Standby/ON-OFF switch illuminates white.

- ▶ Since the connector of the power input socket is used as a disconnect device, the disconnect device must remain freely accessible.

CAUTION!

Protective ground connector!

To avoid the risk of electric shock, this device may only be connected to a power grid which is provided with a protective ground conductor (IEC 60601-1).


- ▶ Only connect the device to power outlets which are equipped with a properly connected protective ground conductor.

NOTE

Risk of tripping!

Inappropriately laid cables represent an increased risk of tripping and falling.

- ▶ Always route cables in such a manner that they do not obstruct your work.

1	"Device Power On/Off" operating button	2	Power input socket
3	Potential equalization connector		

Action

1. Press the yellow locking lever on the device-side connector and slide the connector of the power cord into the power input socket
(exception: The power cord which is used in The Netherlands and Switzerland does not have a locking mechanism).
⇒ When the connector has been fully inserted and you let go of the locking lever, the power cord is secured against unintentional disconnection.
2. Plug the other connector of the power cord into the power outlet.
⇒ The "Device Power On/Off" operating button is illuminated white.
3. If required: connect the potential equalization connector to the equipotential bonding bar in the OR.

4.8 Switching the device on

1	"Device Power On/Off" operating button
---	--

Prerequisite

- ☒ The device is connected to the power supply [124]

Action

1. If you wish to use the following options, before switching the device on, make sure that either the LAN cable is connected or WLAN access to the hospital network has been obtained.
 - ZEISS Smart Services
 - Data storage on a jointly used network drive
 - Data exchange via DICOM
 - Connection of a navigation system via LAN directly or via clinic LAN
2. Push the "Switch system on/off" control key once (max. 1 second).
 - ⇒ The "Device Power On/Off" operating button is illuminated green.
 - ⇒ The operating system is booted.
 - ⇒ The device comes on and displays the start screen and is ready for operation.

* This option is not available in all countries.

4.9 Initial device login


When you use the device for the first time, only the selections "DefaultUser" and "Admin" appear in the "Select User" selection window on the login screen. The "DefaultUser" user group has only limited rights. In order to be able to use the device to a full extent, you must create a new user during initial startup of the device.

Prerequisite

- ☑ The device is connected to the power supply [124] and switched on [126]

Action

1. Tap on the "DefaultUser" field in the "Select User" selection window.

- ⇒ The start screen opens.
- ⇒ In the status bar at the top of the screen, the "DefaultUser" appears in the "User" field.

2. Tap on the [DefaultUser] button.

- ⇒ The "Select User" selection window opens again.
- ⇒ The [Add User] button appears in the "Select User" selection window .

3. Tap on the [Add User] button  and create a new user [201].

- ✓ The next time you start the device, the user you have created will appear in the "Users" selection window. Now you can add new users [201], remove users [204], edit user data [202] etc.

Result

4.10 Connecting the wired foot control panel or rocker foot switch


Figure 72: Connection sockets for foot control panel and rocker foot switch

1	Connection socket for foot control panel	2	Connection socket for rocker foot switch
---	--	---	--

Action

1. If you would like to connect a foot control panel, connect the cable of the foot control panel to the corresponding port on the device.
2. To check the preconfigured button assignment of the foot control panel, open the additional menu "Extended operation" on the right edge of the monitor screen and tap in the [Foot control panel] button .
 - ⇒ The preconfigured button assignment is displayed.
3. To change the button assignment, tap on the [Change settings] button in the display window. Adapt the button assignment as required [181].
4. If you would like to connect a foot control panel, connect the cable of the rocker foot switch to the corresponding port for foot rocker switches on the device.
5. To check the preconfigured button assignment of the rocker foot switch, tap on  Settings →  Operation →  Rocker foot switch. Adapt the button assignment as required [182].

4.11 Pairing the wireless foot control panel with the device


You can use the "Pairing" function to assign the wireless foot control panel (FCB WL) to the device.

Prerequisite

- ☒ Wireless foot control panel is switched on


Action

1. Tap on the [Status Information] button in the status bar.

⇒ The "Status Information" menu is displayed.

2. Tap on the [Pairing] button.

⇒ The "Prepare pairing" menu is displayed.

3. Follow the prompts on the screen. Put the foot control panel in a vertical position in the immediate vicinity of the device.

4. Tap on the [Start] button.

Prepare Pairing

Put the FCP in vertical (rest) position directly next to the device. Press "Start".

⇒ The message "Pairing of the device and the foot control panel in progress" is displayed.

5. "Press a button on the foot control panel until the third LED from the top flashes orange. Keep the foot control panel in a vertical position".

⇒ Following successful pairing, a corresponding message appears on the screen: "Pairing was successfully completed. Place the foot control panel in a horizontal position and perform a function test. Use the wheel on the foot control panel to set the number specified on the stand."

6. Perform a function test. To do this, press any two buttons on the foot control panel.

⇒ The "Radio Connection Intensity" status display lights up green for approx. 1 second.

7. Use the rotary indicator on the foot control panel to set the number specified on the stand sticker and make the assignment of the paired foot control panel to this stand visible.

 CAUTION!

Malfunction due to incorrect pairing!

If pairing is performed incorrectly or not at all, the foot control panel may be disabled or activation of a control may trigger functions on a different stand not assigned to the foot control panel.

- Check and label the assignments of the device and the foot control panel.

8. If pairing is unsuccessful, the "Radio Connection Intensity" status display lights up red for approx. 1s and one of the following messages appears:
"Pairing could not be performed successfully. No foot control panel was detected."
"Pairing could not be performed successfully. Several foot control panels were detected."
9. Repeat the "Pairing" procedure.

4.12 Connecting video systems and external monitors

CAUTION!

Risk of injury caused by electrical voltage.

Connecting unrecognized, defective or non-permitted accessories to live ports may result in electric shock.

- ▶ Connect only recognized, functioning and permitted accessories to live ports.
- ▶ When configuring your ME system, please make sure to comply with the requirements of IEC 60601-1.
- ▶ You must also comply with the Medical Device Operator Ordinance or other national laws.

You can connect the following external monitors to the device:

- Sony 55" 3D* monitor
 - Sony 55" 4K 3D monitor
 - Sony 55" HD 3D monitor
- Sony 55" 2D* monitor
 - Sony 55" 4K 2D monitor
 - Sony 55" HD 2D monitor
- Monitor of a third-party manufacturer [148]

* Sony monitors are not available in all countries.

The video ports are consecutively numbered from 1 to 16 on the device. The following ports are optionally available:

- 1: Available for 4K 2D option
- 1 and 2: Available for 4K 3D option
- 11 to 16: Available for 4K 2D / 3D option


Figure 73: Video ports on the device

1	HDMI video output, left	2	HDMI video output, right
3	DP video input for navigation (MultiVision)	4	DP video input for navigation (MultiVision)
5	DVI video input for external video source**	6	DVI video output for external monitor
7	DVI camera signal video output (touchscreen live image)	8	DVI camera signal video output (touchscreen live image)
9	3G-SDI camera video output	10	3G-SDI camera video output
11	3G-SDI video output, left (2D/3D)	12	3G-SDI video output, left (2D/3D)
13	3G-SDI video output, left (2D)	14	3G-SDI video output, left (2D)
15	3G-SDI video output, right (3D)	16	3G-SDI video output, right (3D)

** Only external video sources with 16:9 format can be used.

4.12.1 Connecting video devices

Approved video cables are listed in product overview G-30-1888.

Action

1. Connect the external video devices to the corresponding video inputs and outputs of the device.
2. Connect the external video sources (e.g. the endoscope) to the corresponding video input of the device.
3. Read and observe the manufacturer's documentation in each case.

4.12.2 Connection of the 4K 3D monitor

You can connect the Sony monitor to the device via the SDI connection sockets or via the HDMI connection sockets.


Figure 74: Connection options for Sony 55" 4K 3D monitor (LMD-X550MT), shown without HDMI connection

Connector ports	Device port	Monitor port	Video cable	Output signals
SDI	11 12 15 16	SDI 1-A SDI 1-B SDI 1-C SDI 1-D	Quad SDI cable, 10 m Order no.: 302584-8785-000	QuadSDI 3840x2160 50/60p (4xBNC)
HDMI	1 or 2	HDMI	HDMI cable, 5 m Order no.: 302584-8758-000	HDMI 2.0 3840x2160 50/60p (ODU)
DVI	6	DVI	DVI cable, 5 m Order no.: 305989-8698-000 DVI cable, 10 m Order no.: 305989-8667-000	DVI-D 1920x1080 50/60p (Lemo)

Table 3: Connection options for Sony 55" 4K 3D monitor

Connecting the device to the Sony 55" 4K 3D monitor

Prerequisite

- ☒ The device and Sony monitor are switched off.
- ☒ An approved video cable is available.

Action

1. Connect the system to the Sony monitor via the SDI connection sockets using 4K-Quad SDI cable.
2. Also connect the device to the Sony monitor via the DVI connection socket.
3. Set the Sony monitor to your desired display mode.

Setting the ZEISS KINEVO 900

4. Switch the device on.
5. Tap on  Settings →  Displays [87].
6. Activate the "Application-specific" function in the "External monitor output" field.
7. Activate the "3D" function in the "Video format" field.

Setting the Sony 55" 4K 3D monitor


Figure 75: Sony 55" 4K 3D monitor control elements

1	"Switch monitor on/off" control key	2	Setting functions
---	-------------------------------------	---	-------------------

8. Push the green main switch on the top of the stand pillar of the monitor cart once.
⇒ The power supply for the Sony monitor is switched on.
9. Push the black "Switch monitor on/off" control key on the lower left edge of the monitor once.
⇒ The Sony monitor is switched on.
10. Push the [Control] button on the lower right edge of the monitor.
⇒ The setting functions are activated.
11. Press the [PORT A] button at the bottom edge of the Sony monitor once or twice.
⇒ Once the [PORT A] button is lit up in green, it is activated.
12. Set the input signal to "SDI1" using the UP and DOWN arrows.
13. Press the [PORT A] button again to exit the menu.
14. Press the [PORT B] button twice.
⇒ Once the [PORT B] button is lit up in green, it is activated.
15. Set the input signal at "DVI-D" using the UP and DOWN arrow keys.
16. Push the [PORT B] button again to exit the menu.
17. Push the [Menu] button once.
⇒ Once the [Menu] button is lit up in green, it is activated.
18. Select the "Screen Control" menu using the UP and DOWN arrows.

19. Push the [Select] button once.

20. Open the "SD11" register using the LEFT and RIGHT arrows and set the following parameters:

SD1 1	SD1 2	DVI-D	HDMI
4K Scan Size		ON	
HD Scan Size		ON	
SD Scan Size		ON	
Zoom		ON	
Flip Pattern		ON	
SD Aspect		4:3	
Interface Mode		4K Quad	

21. Set the "Flip Pattern" function optionally to 180° (Face-to-Face).

22. Open the "DVI-D" register using the LEFT and RIGHT arrows and set the following parameters:

SD1 1	SD1 2	DVI-D	HDMI
HD Scan Size		ON	
SD Scan Size		ON	
Zoom		ON	
Flip Pattern		ON	
SD Aspect		4:3	

23. Push the [Menu] button once.

24. Select the "3D Setting" menu.

25. Open the "DVI-D" tab and set the following parameters:

SD1 1	SD1 2	DVI-D	HDMI
3D/SD Select		3D	
3D Signal Format		Line by Line	
3D Display		0	
L/R Priority		R Line First	

26. Push the [Menu] button again to exit the menu.

Note:

The live image display of the INFRARED 800 with FLOW 800 Option fluorescence application and the QEVO digital exploration tool are superimposed by the white light image of the surgical microscope on the 4K 3D monitor. View these two applications only on the monitors attached to the device.

4.12.3 Connection of the 4K 2D monitor

You can connect the Sony monitor to the device via the SDI connection sockets or via the HDMI connection sockets.


Figure 76: Connection options for Sony 55" 4K 2D monitor (LMD-X550MD), shown without HDMI connection

Connector ports	Device port	Monitor port	Video cable	Output signals
SDI	11 12 13 14	SDI 1-A SDI 1-B SDI 1-C SDI 1-D	Quad SDI cable, 10 m Order no.: 302584-8785-000	QuadSDI 3840x2160 50/60p (4xBNC)
HDMI	1 or 2	HDMI	HDMI cable, 5 m Order no.: 302584-8758-000	HDMI 2.0 3840x2160 50/60p (ODU)
DVI	6	DVI	DVI cable, 5 m Order no.: 305989-8698-000 DVI cable, 10 m Order no.: 305989-8667-000	DVI-D 1920x1080 50/60p (Lemo)

Table 4: Connection options for Sony 55" 4K 2D monitor

Prerequisite

Connecting the device to the Sony 55" 4K 2D monitor

- ☒ The device and Sony monitor are switched off.
- ☒ An approved video cable is available.

Action

1. Connect the device to the Sony monitor either via the SDI connection sockets or via the HDMI connection sockets.
2. Also connect the device to the Sony monitor via the DVI connection socket.

Setting the ZEISS KINEVO 900

3. Switch the device on.
4. Tap on  Settings →  Displays.
5. Activate the "Application-specific" function in the "External monitor output" field.
6. Activate the "2D" function in the "Video format" field.

Setting the Sony 55" 4K 2D monitor


Figure 77: Sony 55" 4K 2D monitor control elements

1	"Switch monitor on/off" control key	2	Setting functions
---	-------------------------------------	---	-------------------

7. Push the green main switch on the top of the stand pillar of the monitor cart once.
 - ⇒ The power supply for the Sony monitor is switched on.
8. Push the black "Switch monitor on/off" control key on the lower left edge of the monitor once.
 - ⇒ The Sony monitor is switched on.

9. Push the [Control] button on the lower right edge of the monitor.
⇒ The setting functions are activated.
10. Press the [PORT A] button at the bottom edge of the Sony monitor once or twice.
⇒ Once the [PORT A] button is lit up in green, it is activated.
11. Set the input signal to "SDI1" using the UP and DOWN arrows.
12. Press the [PORT A] button again to exit the menu.
13. Press the [PORT B] button twice.
⇒ Once the [PORT B] button is lit up in green, it is activated.
14. Set the input signal at "DVI-D" using the UP and DOWN arrow keys.
15. Push the [PORT B] button again to exit the menu.
16. Push the [Menu] button once.
⇒ Once the [Menu] button is lit up in green, it is activated.
17. Select the "Screen Control" menu using the UP and DOWN arrows.
18. Push the [Select] button once.
19. Open the "SDI1" register using the LEFT and RIGHT arrows and set the following parameters:

SDI 1	SDI 2	DVI-D	HDMI
4K Scan Size		Off	
HD Scan Size		Off	
SD Scan Size		Off	
Zoom		Off	
Flip Pattern		Off	
SD Aspect		4:3	
Interface Mode		4K Quad	


20. Set the "Flip Pattern" function optionally to 180° (Face-to-Face).
21. Open the "DVI-D" register using the LEFT and RIGHT arrows and set the following parameters:

SDI 1	SDI 2	DVI-D	HDMI
HD Scan Size		Off	
SD Scan Size		Off	
Zoom		Off	
Flip Pattern		Off	
SD Aspect		4:3	

22. Push the [Menu] button once.
23. Select the "3D Setting" menu.
24. Open the "SDI1" tab and set the following parameters:

SDI 1	SDI 2	DVI-D	HDMI
3D/2D Select		2D	
3D Signal Format			
3D Clarity			
L/R Priority			

25. Open the "DVI-D" tab and set the following parameters:

26. Push the [Menu] button again to exit the menu.

Note:

The live image display of the INFRARED 800 with FLOW 800 Option fluorescence application and the QEVO digital exploration tool are superimposed by the white light image of the surgical microscope on the 4K 3D monitor. View these two applications only on the monitors attached to the device.

4.12.4 Connection of the HD 3D monitor


Figure 78: Connection options for Sony 55" HD 3D monitor

Connector ports	Device port	Monitor port	Video cable	Output signals
SDI	9 10	SDI 1-A SDI 1-B	SDI cable, 10 m Order no.: 305989-8764-000	3G-SDI or HD-SDI 1920x1080 50/60i/p ("Interlace" or "Progressive", selectable), 2xBNC
DVI	6	DVI	DVI cable, 5 m Order no.: 305989-8698-000 DVI cable, 10 m Order no.: 305989-8667-000	DVI-D 1920x1080 50/60p (Lemo)

Table 5: Connection options for Sony 55" HD 3D monitor

Connecting the device to the Sony 55" HD 3D monitor

Prerequisite

- ☒ The device and Sony monitor are switched off.
- ☒ An approved video cable is available.

Action

1. Connect the device to the Sony monitor via the SDI connection sockets.
2. Also connect the device to the Sony monitor via the DVI connection socket.

Setting the ZEISS KINEVO 900

3. Switch the device on.
4. Tap on  Settings →  Displays [87].
5. Activate the "Application-specific" function in the "External monitor output" field.
6. Activate the "3D" function in the "Video format" field.

Setting the Sony 55" HD 3D monitor


Figure 79: Sony Monitor 55" HD 3D control elements

1	"Switch monitor on/off" control key	2	Setting functions
---	-------------------------------------	---	-------------------

7. Push the green main switch on the top of the stand pillar of the monitor cart once.
⇒ The power supply for the Sony monitor is switched on.
8. Push the black "Switch monitor on/off" control key on the lower left edge of the monitor once.
⇒ The Sony monitor is switched on.
9. Push the [Control] button on the lower right edge of the monitor.
⇒ The setting functions are activated.
10. Press the [PORT A] button at the bottom edge of the Sony monitor once or twice.
⇒ Once the [PORT A] button is lit up in green, it is activated.
11. Set the input signal to "SDI1" using the UP and DOWN arrows.
12. Press the [PORT A] button again to exit the menu.
13. Press the [PORT B] button twice.
⇒ Once the [PORT B] button is lit up in green, it is activated.

14. Set the input signal at "DVI-D" using the UP and DOWN arrow keys.
15. Push the [PORT B] button again to exit the menu.
16. Push the [Menu] button once.
⇒ Once the [Menu] button is lit up in green, it is activated.
17. Select the "Screen Control" menu using the UP and DOWN arrows.
18. Push the [Select] button once.
19. Open the "SDI1" register using the LEFT and RIGHT arrows and set the following parameters:

SDI 1	SDI 2	DVI-D	HDMI
4K Scan Size		Off	
HD Scan Size		Off	
SD Scan Size		Off	
Zoom		Off	
Flip Pattern		Off	
SD Aspect		4:3	
Interface Mode		HD Dual	

20. Set the "Flip Pattern" function optionally to 180° (Face-to-Face).
21. Open the "DVI-D" register using the LEFT and RIGHT arrows and set the following parameters:

SDI 1	SDI 2	DVI-D	HDMI
HD Scan Size		Off	
SD Scan Size		Off	
Zoom		Off	
Flip Pattern		Off	
SD Aspect		4:3	

22. Push the [Menu] button once.
23. Select the "3D Setting" menu.
24. Open the "SDI1" tab and set the following parameters:

SDI 1	SDI 2	DVI-D	HDMI
2D/3D Select		3D	
3D Signal Format		Dual Stream	
3D Disparity		0	
L/R Priority		L Line First	

25. Open the "DVI-D" tab and set the following parameters:

SDI 1	SDI 2	DVI-D	HDMI
2D/3D Select		2D	
3D Signal Format		Line by Line	
3D Disparity		0	
L/R Priority		R Line First	

26. Push the [Menu] button again to exit the menu.

Note:

The live image display of the INFRARED 800 with FLOW 800 Option fluorescence application and the QEVO digital exploration tool are superimposed by the white light image of the surgical microscope on the 4K 3D monitor. View these two applications only on the monitors attached to the device.

4.12.5 Connection of the HD 2D monitor


Figure 80: Connection options for Sony 55" HD 2D monitor

Connector ports	Device port	Monitor port	Video cable	Output signals
SDI	9 or 10	SDI 1-A	SDI cable, 10 m Order no.: 305989-8764-000	3G-SDI or HD-SDI 1920x1080 50/60i/p ("Interlace" or "Progressive", selectable), 2xBNC
DVI	6	DVI	DVI cable, 5 m Order no.: 305989-8698-000 DVI cable, 10 m Order no.: 305989-8667-000	DVI-D 1920x1080 50/60p (Lemo)

Table 6: Connection options for Sony 55" HD 2D monitor

Connecting the device to the Sony 55" HD 2D monitor

Prerequisite

- ☑ The device and Sony monitor are switched off.
- ☑ An approved video cable is available.

Action

1. Also connect the device to the Sony monitor via the SDI connection socket.
2. Also connect the device to the Sony monitor via the DVI connection socket.

Setting the ZEISS KINEVO 900

- 3. Switch the device on.
- 4. Tap on  Settings →  Displays.
- 5. Activate the "Application-specific" function in the "External monitor output" field.
- 6. Activate the "2D" function in the "Video format" field.

Setting the Sony 55" HD 2D monitor


Figure 81: Sony Monitor 55" HD 2D control elements

1	"Switch monitor on/off" control key	2	Setting functions
---	-------------------------------------	---	-------------------

- 7. Push the green main switch on the top of the stand pillar of the monitor cart once.
⇒ The power supply for the Sony monitor is switched on.
- 8. Push the black "Switch monitor on/off" control key on the lower left edge of the monitor once.
⇒ The Sony monitor is switched on.
- 9. Push the [Control] button on the lower right edge of the monitor.
⇒ The setting functions are activated.
- 10. Press the [PORT A] button at the bottom edge of the Sony monitor once or twice.
⇒ Once the [PORT A] button is lit up in green, it is activated.
- 11. Set the input signal to "SDI1" using the UP and DOWN arrows.
- 12. Press the [PORT A] button again to exit the menu.
- 13. Press the [PORT B] button twice.
⇒ Once the [PORT B] button is lit up in green, it is activated.

14. Set the input signal at "DVI-D" using the UP and DOWN arrow keys.
15. Push the [PORT B] button again to exit the menu.
16. Push the [Menu] button once.
⇒ Once the [Menu] button is lit up in green, it is activated.
17. Select the "Screen Control" menu using the UP and DOWN arrows.
18. Push the [Select] button once.
19. Open the "SDI1" register using the LEFT and RIGHT arrows and set the following parameters:

SDI 1	SDI 2	DVI-D	HDMI
4K Scan Size		Off	
HD Scan Size		Off	
SD Scan Size		Off	
Zoom		Off	
Flip Pattern		Off	
SD Aspect		4:3	
Interface Mode		HD/SD Single	

20. Set the "Flip Pattern" function optionally to 180° (Face-to-Face).
21. Open the "DVI-D" register using the LEFT and RIGHT arrows and set the following parameters:

SDI 1	SDI 2	DVI-D	HDMI
HD Scan Size		Off	
SD Scan Size		Off	
Zoom		Off	
Flip Pattern		Off	
SD Aspect		4:3	

22. Push the [Menu] button once.
23. Select the "3D Setting" menu.
24. Open the "SDI1" tab and set the following parameters:

SDI 1	SDI 2	DVI-D	HDMI
2D/3D Select		2D	
3D Signal Format			
3D Disparity			
L/R Priority			

25. Open the "DVI-D" tab and set the following parameters:

SDI 1	SDI 2	DVI-D	HDMI
2D/3D Select		2D	
3D Signal Format			
3D Disparity			
L/R Priority			

26. Push the [Menu] button again to exit the menu.

Note:

The live image display of the INFRARED 800 with FLOW 800 Option fluorescence application and the QEVO digital exploration tool are superimposed by the white light image of the surgical microscope on the 4K 3D monitor. View these two applications only on the monitors attached to the device.

4.12.6 Connecting the monitor of a third-party manufacturer

Prerequisite

- ☒ The device and third-party-manufacturer monitor are switched off

Action

1. Read and observe the monitor manufacturer's documentation.
2. Connect the monitor to the corresponding connection sockets of the device [47] with approved video cables.
3. Set the monitor to your desired display mode. Read and observe the monitor manufacturer's documentation for this purpose.

4.13 Connecting USB storage media

NOTE

Damage resulting from connecting the device to an external storage device

Connecting the device to an external storage device may permit viruses to be transferred to the device or may grant viruses unauthorized access to the system.

- ▶ Use only USB storage devices which have been approved by ZEISS.
- ▶ Make sure that the USB storage device is free of viruses. The operator is responsible for security.
- ▶ Protect the network by using suitable safeguards (e.g. firewalls) against unauthorized access.

Connecting a USB storage device

Note:

The device can detect only one USB storage device. For this reason, only connect the storage device that you want to save your data to.

Prerequisite

- ☒ All USB storage devices which are not currently needed have been removed.

Action

1. If applicable, press the push button to open the USB cover [46].
2. Connect the USB storage device to one of the device's USB ports.
 - ⇒ The storage device will be recognized automatically.
3. Place the storage device in the provided tray if needed.

Removing a USB storage device

Note:

To prevent your data from being lost, you must eject the USB storage device before removing it.

4. Tap on the [Status Information] button in the status bar.
⇒ The "Status Information" menu is displayed.
5. In the "Storage" section under "USB", tap [Eject].
⇒ The storage device will be safely disconnected from the device.
6. Remove the USB storage device.
7. If applicable, close the USB cover.

4.14 Configure network

NOTE

Damage resulting from connecting the device to an external storage device (USB storage device) or an external network.

Connecting the device to external storage devices or networks may permit viruses to be transferred to the device or may grant viruses unauthorized access to the system.

- ▶ Use only ZEISS-approved USB storage media.
- ▶ Make sure that the USB storage device is free of viruses. The operator is responsible for security.
- ▶ Protect the network by using suitable safeguards (e.g., firewalls) against unauthorized access.

You can connect the device to your hospital network in order to export patient data to your network or to a DICOM server or in order to import patient data from a DIOCOM worklist.

The following settings and configurations may only be carried out by authorized users with IT admin system rights [149].

4.14.1 Activating IT admin system rights

The IT admin password is contained in a sealed red envelope included in the scope of supply.

If you have forgotten the password, contact ZEISS Service. You can find the ZEISS contact person for your country on the Internet at the following website: www.zeiss.com/med.


Action

1. Tap on the [User] button in the status bar.
⇒ The "Select User" menu is displayed. The current user is highlighted in blue.
2. Select the "Admin" user.
⇒ A window for entering the password opens.

3. Enter the IT admin password.
4. Tap on the [OK] button .
⇒ The "Admin" user is displayed in the status bar.

4.14.2 Activate network connection via LAN

The Ethernet port (LAN) on the device contains a network isolator as per IEC 60601-1 that galvanically decouples the Ethernet connection between the network and the device.

1	Network/Ethernet connection (LAN)
---	-----------------------------------

NOTE

Conflicting addresses

Some IP addresses areas are reserved for certain functionalities and may not be used to establish a network.

- ▶ Do not enter IP addresses from the **192.168.173.xxx** area (reserved for hotspot).
- ▶ Do not enter an IP address from the **192.0.2.xxx** area (reserved for QEVO connection (Option)).

Prerequisite

- ☒ The device is connected to the respective network via a network cable.

Action

1. Tap on ⚙ Settings → Connections.
2. Scroll to the paragraph titled "LAN".
3. If necessary, enter network data.
4. Tap the [Save] button.
⇒ The network data is stored.

5. Activate the "Connect with LAN" button.

Result

- ✓ The LAN connection is active.

4.14.3 Activate network connection via wireless LAN (option)

The WLAN option is not available in all countries.

NOTE

Conflicting addresses

Some IP addresses areas are reserved for certain functionalities and may not be used to establish a network.

- ▶ Do not enter IP addresses from the **192.168.173.xxx** area (reserved for hotspot).
- ▶ Do not enter an IP address from the **192.0.2.xxx** area (reserved for QEVO connection (Option)).

Prerequisite

- ☒ The network corresponds to the WPA2/PSK or WPA2/PEAP standards.
- ☒ If necessary: The PEAP authentication certificate to be loaded has the ending ".cer", ".der" or ".pem" and is in the root directory of a connected USB storage medium.

Action

1. Tap on  Settings → Connections.
2. Scroll to the paragraph titled "WLAN".
3. In the "Authentication" selection field, select the desired authentication standard.
4. Tap the [Import certificate] button to load a PEAP authentication certificate.
5. If necessary, insert a check mark in the "DHCP" field or enter the network data.
6. Tap the [Save] button.
 - ⇒ The network data is stored.
7. Activate the "Connect with WLAN" button.

Result

- ✓ The WLAN connection is active.

4.14.4 Creating and activating the device as a hotspot

The device has to be created and activated as a hotspot so that the "Default Users" and "Users" user groups can use the device as a hotspot.

The WLAN option required to use the device as a hotspot is not available in all countries.

Action

1. Tap on  Settings → Connections.
2. Scroll to the paragraph titled "Hotspot".
3. If necessary, enter the name and password of the hotspot.
4. Tap the [Save] button.
 - ⇒ The hotspot data is saved.

5. Insert the checkmark for "Activate hotspot".

Result

- ✓ The device is activated as a hotspot and can be used as a hotspot [239] by the "Default Users" and "Users" user groups.

4.14.5 Activating the web interface in the device

If the device is created and activated as a hotspot [151], the patient data may be accessed via an external device. A password must be allocated by the IT administrator for this purpose.

The WLAN option required to use the device as a hotspot is not available in all countries.

Prerequisite

- ☒ The device is created as a hotspot and activated [151].

Action

1. Tap on  Settings → Connections.
2. Scroll to the paragraph titled "Web interfaces".
3. Enter the password for the web interface.
4. Tap the [Save] button.

Result

- ✓ The "Default Users" or "Users" user group can access patient data using an external device [240] via the hotspot after entering the password.

4.14.6 Network drive

The device enables the storage of patient data (images and videos) in a shared network drive folder.

Prerequisite

- ☒ The hospital server supports the "Server Message Block (SMB)" and/or "Common Internet File System (CIFS)" protocols.
- ☒ The device is connected to the respective network via LAN [150] or WLAN [151].

Action

1. Tap on  Settings → Connections.
2. Scroll to the paragraph titled "Network drive".
3. If necessary, enter the corresponding network drive data.
4. Tap the [Save] button.
 - ⇒ The network drive data is saved.
5. Activate the "Connect with network approval" button.

Result

- ✓ The connection to the network drive is active.

4.14.7 Configure the network connection for DICOM (option).

Network and DICOM addresses must be assigned by the IT administrator in charge. ZEISS Service must be aware of all IT and DICOM parameters (station name, AE title, IP address, port number) ahead of the installation.

NOTE

Incorrect or incomplete configuration

Inadequate knowledge of the DICOM configuration requirements may lead to defective settings and connection errors.

- During network setting, read and observe the ZEISS DICOM Conformance Statement (G-30-1952).

Prerequisite

- ☒ The device is connected to the network via LAN [150] or WLAN [151].

Action

1. Tap on  Settings → Connections.
2. Scroll to the paragraph titled "DICOM".
3. Enter all DICOM data. Preset data may be overwritten as necessary.
NOTE! Please pay attention to upper-case and lower-case letters.
4. Tap the [Save] button.
5. Tap on the [Test Connection] button.
⇒ The result of the connection test appears in an info window.
6. If the connection test was unsuccessful, check the settings made and the network connection.

You can make other settings to your DICOM connection in the "Advanced DICOM settings" paragraph, as necessary.

Prerequisite

- ☒ The Settings → Connections menu is open.

7. Scroll to the paragraph titled "Advanced DICOM settings".
8. Enter the desired data and select the desired options.
9. Tap the [Save] button.
⇒ The settings are saved.

4.14.8 Changing the computer name of the device

Action

1. Tap on  Settings → Connections.
2. Scroll to the paragraph titled "Computer name".
3. Enter the computer name in the "Change PC name" field.
4. Tap the [Save] button.
⇒ A system notification appears on the monitor.
5. Confirm the system notification.
6. Reboot the device.

4.15 Configuring the Service PC

The following settings and configurations may only be carried out by authorized users with IT admin system rights [149].

You can adjust the following settings:

- **Service**
 - Synchronize Patient Data
 - Export User
 - Import User
- **Settings**
 - Activate User Password
 - Change standard password
 - Set System Language
 - Change video frequency
 - Setting the date and time
 - Change Time Format
 - Change Date Format
- **Log Files**
 - Export Log Files
 - Export Audit Files
 - Load DICOM Log Files

4.15.1 Service

Action

1. Tap on ⚙ Settings → Service PC → Service.
2. Tap on the corresponding button [Start] in order to initiate the action.
3. Close the "Settings" menu by tapping on the ↗ button.

4.15.2 Settings

4.15.2.1 Activating and changing the user password

If you enable the user password function in order to protect patient data, all users must log onto the device using the standard password "zeisszeiss". All newly created users are also assigned the standard password. The standard password can be changed by the IT admin. The IT admin is responsible for the standard password and for providing it to all users of the device.

Users can log in to the device as a "Default User"; no password needs to be entered to do this. However, they will then only be able to view the current patient's data.

Action

1. Tap on  Settings → Service PC → Settings.
2. Tap on the slide switch in the "Activate User Password" field to activate the user password.
 - ⇒ The slide switch will turn blue .
 - ⇒ The user password is activated.
 - ⇒ All users must log in on the device using the standard password that has been issued to them.
 - ⇒ Newly added users are automatically assigned the issued standard password.
3. In the "Standard Password" field, tap on the [Change] button to change the standard password.
4. Enter a new password.
5. Confirm your entry by tapping on the [Done] button.
6. Disclose the new user password to the users of the device.
 - ⇒ All users must log in on the device using the standard password that has been issued to them.
 - ⇒ Newly added users are automatically assigned the issued standard password.

Note:

Disabling and re-enabling the user password will cause the passwords of all users to be reset to the standard password.

4.15.2.2 Setting general parameters

Action

1. Tap on  Settings → Service PC → Settings.

You can set the following general parameters:

2. **Set System Language:** → Select the desired system language.
3. **Change Video Frequency:** → 50 Hz or 60 Hz.
 - ⇒ The changes take effect only after the device is restarted.
4. **Set Date and Time:** → After setting, tap on the [Accept] button.
5. **Change Time Format:** → 24h or 12h.
6. **Change Date Format:** → DD.MM.YYYY or MM/DD/YYYY.
7. Close the "Settings" menu by tapping on the  button.

4.15.3 Log Files

Action

1. Tap on  Settings → Service PC → Log Files.
2. Tap on the corresponding button [Start] in order to initiate the action.
3. Tap on the button for the file type and select the file type to be loaded.
4. Tap on the [Load log file] button to start the action.
5. Close the "Settings" menu by tapping on the  button.

4.16 Connecting the navigation system to the device

You can connect navigation systems to the device:

- Navigation system via standard connection (11). The motorized positioning and travel functions of the navigation system are limited to the three microscope axes.
- Navigation system via extended connection (network LAN). The motorized positioning and travel functions of the navigation system are possible in all stand and microscope axes.

The initial connection (calibration) of an external navigation system must be performed by trained personnel of the navigation system manufacturer. Please observe the operating instructions provided in the Instructions for Use of the navigation system.

CAUTION!

Faulty navigation!

The navigation antenna attached to the microscope may become maladjusted due to a collision during transport or coarse positioning.

- ▶ Check the entire navigation system according to the manufacturer's specifications before every use.
- ▶ Check the calibration of the navigation system according to the manufacturer's specifications before every use.
- ▶ Before every intervention using a connected and authenticated navigation system, the function and accuracy of the navigation system including the display in the data injection must be verified accordingly (e.g. by focusing on a measuring point or comparing the focal point with a navigated instrument). Observe the corresponding specifications of the navigation system.
- ▶ Connection of the navigation system to the device results in a medical system for which the system supplier (manufacturer of the navigation system) must meet the stipulated requirements (approval, qualifications, etc.). All accompanying papers required will be supplied by the manufacturer of the navigation system.

Prerequisite

- ☑ The navigation antenna has been installed on the microscope (by trained personnel from the navigation system manufacturer)

Action

1. Connect the navigation system to the device; for this purpose, read and observe the navigation system manufacturer's corresponding instructions.
2. To use the data injection functions of the navigation system, activate the navigation system in the "Settings MultiVision" menu [168].

4.17 Mounting the laser micromanipulators

CAUTION!

Incorrect mounting of a laser micromanipulator

Incorrect mounting of the laser micromanipulator can cause damage to equipment and errors during use. Anyone connecting additional equipment to medical electrical systems is a system configurator and as such responsible for compliance of the system with the standards for systems (see IEC 60601-1). Local legislation has priority over the requirements of the abovementioned standards.

- Observe the information and the instructions for use supplied by the manufacturer of the laser micromanipulator. If necessary, contact your ZEISS representative.
- Before fitting the ZEISS dovetail mount, consult the laser manufacturer to check whether it is compatible. Some laser manufacturers offer their own assembly sets for their laser micromanipulators.


Figure 82: Fitting the dovetail mount on the underside of the surgical microscope

1	Threaded bores (4x) for fastening the dovetail mount for laser micromanipulators
---	--

Prerequisite

- ☒ The microscope is configured as a system for hybrid (optical-digital) visualization.
- ☒ Dovetail mount including four M3x12 fastening screws (screw material: A2-70) is available (order number: 303360-9903-000).
- ☒ The laser micro manipulator weighs no more than 1 kg.
- ☒ The system is switched off.

Action

1. Fasten the dovetail mount for laser micromanipulators to the bottom of the surgical microscope using the four (M3x12) screws supplied.
2. Make sure that the maximum permissible screw-in depth of the fastening screws into the Varioskop carrier is 5 mm.
3. Tighten the four screws firmly to a tightening torque of 0.93 Nm.
4. Slide the laser micromanipulator into the dovetail mount.
5. Tighten the knurled screw firmly to secure the laser micromanipulator in the dovetail mount.

Result

- ✓ You can prepare the device for the use of a laser micromanipulator.

5 Daily startup

5.1 Safety during preparation

CAUTION!

Functional deterioration!

We recommend taking adequate precautions, depending on the application, to enable the surgical procedure or treatment to be finished without using this microscope (for example in case of a system defect).

As the age of the light source increases, the illumination intensity achieved in a particular setting decreases. Please replace the xenon lamp in due time. Please note the remaining service hours displayed on the touchscreen.

- ▶ Check the instrument for proper functionality or damage before each use.
- ▶ Perform a function test; if required, replace the lamps prior to using the system.
- ▶ Make sure that no tissue damage is caused by excessive illumination intensity.
- ▶ Always keep a backup lamp ready so that you can replace a defective lamp following surgery.

NOTE

Do not cover the ventilation openings!

Overheating of the device

- ▶ Do not cover the ventilation openings, as overheating of the light source will cause the device to switch off.

5.2 Setting the position of the handgrips


Figure 83: Setting the position of the hand grips

1	Locking lever	2	Handle
3	Handgrip holder	4	Handgrip holder locking screw

Action

1. Move the surgical microscope so that you have a good view of the surgical field.
2. Open both locking screws to adapt the position of the handgrip holders.
3. Move the handgrip holders with the handgrips to a position that is favorable and comfortable for you.
4. Tighten both locking screws firmly.
5. If necessary, open the two clamping levers and perform a fine adjustment of the position of both handgrips.
⇒ The handgrips can be rotated approx. 180°.
6. Swivel or rotate the handgrips such that you are able to operate the surgical microscope comfortably.
7. Press the clamping levers shut hand-tight.
8. Then check whether the required swiveling ranges of the surgical microscope are sufficient and any collisions with attached components (e.g. the lateral co-observer) can be excluded.

5.3 Positioning the device in the OR


Figure 84: Positioning the device in the OR

1	Possible positions for cranial procedures	2	Possible positions for face-to-face (spine) procedures
3	Surgeon	4	Assistant
5	Sterile specialist	6	Non-sterile specialist

5.4 Check that the system can move freely!

⚠ CAUTION!

Check that the system can move freely!

Position the device in such a way that you can move it away from the patient at any time by overpowering the brake mechanism. Since the connector of the power input socket is used as a disconnect device, the disconnect device must remain freely accessible.

5.5 Configuring the device for applications

The information in this paragraph only applies to the optional configuration of the microscope as a system for hybrid (optical-digital) visualization.


Figure 85: Configuring the device for applications

1	Configuring the device for "face-to-face" co-observation (spine)	2	Configuring the device for "left/right" co-observation (cranial)
---	--	---	--

5.6 Drape position

The drape position is saved as a factory setting, however, can be changed arbitrarily.

⚠ CAUTION!

Risk of crushing!

Fingers may be crushed between the vertical arm and the horizontal arm.

- ▶ The drape position may be used only without a patient!
- ▶ Never touch the area between the vertical arm and the horizontal arm while you are moving to the drape position.

Action

1. Tap on ⚙ Settings → 🏠 Stand.
2. Tap on the [Start] button in the "Start Drape Position" field.
3. To activate the travel, press the XY joystick button on the right handgrip to the right  and hold it down until the drape position has been reached.
The joystick on the FCP also can be used to activate the travel.
 - ⇒ The device then automatically travels to the drape position.
 - ⇒ When a beep sounds, the drape position has been reached.

Determining a new drape position

If the factory set drape position is unfavorable for you, you can define and save your own drape position.

4. Move the device to a position that is favorable and comfortable for you.
5. Tap on  Settings →  Stand.
6. Tap on the [OK] button in the "Save Current Position as Drape Position" field.
 - ⇒ The current position of the device is saved as the new drape position.
7. Now you can move the device to the drape position as usual.

5.7 Attaching a ZEISS SMARTDRAPE

CAUTION!

Risk of infection

Risk of infection due to contamination.

- For sterile draping of the device, use disposable sterile covers (ZEISS SMARTDRAPE).

CAUTION!

Using drapes from other manufacturers may result in deviations of the optical parameters.

The cover glasses on drapes from other manufacturers may impair the optical characteristics of the device. This may result e.g. in the transmission of falsified navigation data.

- Use only a ZEISS SMARTDRAPE for the device.


Figure 86: Attaching a SMARTDRAPE

1	Non-sterile area	2	Sterile area
3	Effective area of drape vacuum system	4	Marking: Close SMARTDRAPE air-tight here
5	RFID reader		

Prerequisite

- ☑ The device is in the drape position [162]
- ☑ When using drapes from third-party manufacturers without RFID detection, the "AutoDrape" function (air extraction for easier draping) is not available. The magnetic fastening function cannot be used if the drape does not have a corresponding magnetic holder.
- ☑ ZEISS SMARTDRAPEs (order number 306028-0000-000) are tested for use with navigation systems.

Action

1. Unpack the SMARTDRAPE.
2. Pull the SMARTDRAPE over the surgical microscope.
3. Fasten the drape ring of the SMARTDRAPE to the lens below the surgical microscope.
 - ⇒ The drape ring with the magnet engages in the correct position (proper orientation of the beveled cover glass to prevent reflections).
4. Pull the SMARTDRAPE over the horizontal arm until it covers the marking.

When attaching the SMARTDRAPE, make sure that there is sufficient free space for swiveling, tilting and rotating movements of the surgical microscope.

Make sure the control elements can be safely operated through the SMARTDRAPE.

5. Seal the SMARTDRAPE airtight at the marking using one of the fastening straps.
6. Remove the RFID label from the SMARTDRAPE.
7. Stick the adhesive label on the RFID reader.
 - ⇒ It can take up to 10 seconds to successfully read in the code.
 - ⇒ A beep sounds as soon as the code has been successfully read in.
 - ⇒ The AutoDrape function is available.
8. Open the additional menu "Extended operation" on the touch-screen.
9. Open the additional menu "Extended operation" on the touch-screen.
10. Activate the [AutoDrape] slide switch to switch on the automatic suction pump.
 - ⇒ The vacuum system operates at maximum power after approx. 2 minutes.
 - ⇒ After approx. 2 minutes, the suction system automatically switches to a reduced power to maintain a vacuum.
11. Remove the RFID label once the operation has ended.

5.8 Autobalance

Figure 87: Autobalance

CAUTION!

Use the device in perfectly balanced condition only!

With an incorrectly balanced device, actuation of the brake buttons may lead to uncontrolled movements of the device.

- ▶ Balancing and the subsequent test may not be carried out when the system is positioned over the patient and may be done only at a safe distance from other persons and devices.
- ▶ Check the balance of the device by actuating the brake buttons while holding the microscope securely by both hand grips. Repeat the autobalance if necessary.

Prerequisite

- ☒ No tubes and accessories are mounted on the system in the fully digital microscope configuration.

Action

1. Autobalance full system: Tap the [System] button in the "Extended operation" menu on the touchscreen.
2. Autobalance drape (balance out axes 2 and 3): Tap on the [Drape] button in the "Extended operation" additional menu on the touchscreen.

If the device has been correctly balanced, the surgical microscope can be moved almost effortlessly.

5.9 Setting the microscope

The information in this paragraph only applies to the optional configuration of the microscope as a system for hybrid (optical-digital) visualization.

Action

1. Move the microscope to a vertical position above a flat object, e.g. a sheet of paper with writing.
2. Adjust the distance of the eyepieces at the binocular tube to your interocular distance so that the two eyepiece images (object and edge of visual field) merge into one.
3. Set the microscope to minimum magnification. Select a working distance which is as short as possible (shortest working distance + approx. 25 mm).
4. Set the diopter setting ring on the eyepiece to 0 D. (Diopters).
5. Look through the eyepiece and focus the image sharply.
6. Set the microscope to maximum magnification and adjust the fine focusing system until the image is sharply defined.
7. Then reset the microscope to minimum magnification without changing the working distance.
8. Adjust the diopter setting ring on the eyepiece to the maximum positive dpt. value (e.g. +5 dpt.).
9. Look through the eyepiece and slowly turn the diopter setting ring in the minus dpt. direction until the image is once again defined sharply.

10. Repeat the entire procedure for the second eyepiece.

⇒ The microscope has now been set in such a way that the image is sharply defined over the complete magnification range without any renewed focusing during the magnification setting.

11. If you nevertheless need to refocus, repeat this procedure.

12. Adjust the eyecups of the eyepieces in such a way that the entire visual field can be seen.

Observation using eyeglasses: Screw the eyecups in all the way.

Observation without eyeglasses: Unscrew the eyecups and adjust them individually.

5.10 Connect QEVO Micro-Inspection Tool (option)

When connecting and using the Micro-Inspection Tool, it is vital that you observe the relevant safety instructions and warning notices in the QEVO and QEVO ECU Instructions for Use! [13]

Prerequisite

☑ The monitoring and control unit QEVO ECU is installed on ZEISS KINEVO 900 .

☑ QEVO is prepared in sterile form.

☑ QEVO is placed in the sterile reprocessing basketInstru-Safe.

Action

1. Sterile person: Remove QEVO from the reprocessing basket with the cableInstru-Safe.

2. Sterile person: Deposit the QEVO in the sterile area safely.

3. Unsterile person: Plug the connector into the connection socket of the control unit.

4. Lay the cable so that it does not disturb any user and nobody can stumble over it.

In order to use the Micro-Inspection Tool, you must also activate QEVO in the user surface of ZEISS KINEVO 900 . [168]

5.11 Activate QEVO Micro-Inspection Tool (option)

Prerequisite

- ☒ The license for the Micro-Inspection Tool QEVO has been activated.

Action

- ☒ ZEISS KINEVO 900 is switched on.

1. Tap on  Settings →  MultiVision.
2. To enable the MultiVision Data Injection of the QEVO camera image, activate the "QEVO / Video" function in the "Source" field.
 - ⇒ When pressing the [MultiVision] button on the hand grip or the foot control panel [180] the camera image is injected by QEVO.
3. To enable the display of the QEVO camera image on the monitor, activate the "QEVO" function in the "Video source" field ("QEVO" section).
 - ⇒ When pressing the [QEVO] button on the hand grip or foot control panel [180] the QEVO camera image is shown on the system monitors as a full image. The ZEISS KINEVO 900 camera image is shown on the touchscreen as Picture in Picture.

5.12 Activating a connection to the navigation system

Please note the information provided by the manufacturer of the navigation system. The navigation system to be connected must be calibrated once for the ZEISS KINEVO 900 by the manufacturer in order to be used properly.

CAUTION!

Faulty navigation!

The navigation antenna may become maladjusted due to a collision during transport or coarse positioning.

- Check the calibration of the navigation system according to the manufacturer's specifications before every use.

Prerequisite

- ☒ The navigation license has been activated.
- ☒ The navigation system is connected to the device [156]

Action

1. Tap on  Settings →  MultiVision.
2. Activate the "MultiVision" slide switch.
 - ⇒ The "MultiVision" function is switched on, the data injection of the navigation system can be activated.

Result

3. Activate the "Navigation" function in the "Source" field.
4. Depending on the navigation system connected, activate the "Navigation (standard)" or the "Navigation ("extended")" function.
 - ⇒ As soon as you activate the slide switch in the "Navigation (extended)" field, the current IP address appears under the slide switch.
 - ✓ Following successful activation and inspection, you can install the functionalities of the navigation system on the device and use them there.

5.13 Preparing the device for the use of laser micromanipulators

CAUTION!

Danger due to unintentional change to the focus position.

Due to deactivation of the "Focus Lock" function, the focal planes of the surgical microscope and the laser micromanipulator do not match and the desired effectiveness of the laser micromanipulator is not achieved.

- ▶ Before each use, check that the focal planes correspond (without a patient present). It is essential that the focal plane of the surgical microscope matches that of the laser micromanipulator.
- ▶ When operating an approved external laser micromanipulator, activate the "Focus Lock" function if it does not have an AF function.
- ▶ Observe the information and the instructions for use supplied by the manufacturer of the laser micromanipulator. If necessary, contact your ZEISS representative.

Prerequisite

- ☒ The microscope is configured as a system for hybrid (optical-digital) visualization.
- ☒ The laser micromanipulator is mounted on the surgical microscope
- ☒ No patient is lying underneath the surgical microscope and the mounted laser micromanipulator

Action

1. Switch the device on.
2. Adjust the working distance (rough focus) to the focus value of the laser micromanipulator.
 - ⇒ The current focus value is displayed on the monitor.
3. If necessary, correct the focus adjustment by slightly shifting it (fine focus).
4. Aim the aiming beam at a suitable object (e.g. a wooden spatula) to check the mutual alignment of the focal planes.

5. Trigger a single pulse.

⇒ The center of the burn zone must not deviate from the center of the aiming beam by more than max. 0.5 mm.

6. If necessary, perform minor adjustments to correct the focus.

7. Tap on ⚙ Settings → 🔍 Microscope → Focus.

8. Activate the slide switch in the "Focus Lock" field.

⇒ All motorized focus functions, including the "Autofocus" function, are deactivated.

⇒ The focal planes cannot be adjusted by motor.

⇒ The "Autofocus" function will not be triggered when the brakes are activated.

9. In order to work with the laser micromanipulator, read and observe the manufacturer's documentation.

5.14 Device function tests

Use the following list to check device functionality before each use (patient should not be present!).

Action

- ▶ Check AutoBalance [165].
 - ⇒ The surgical microscope maintains its balance in all positions of the working range when the magnetic brakes are released.
- ▶ Check accessories on the surgical microscope [116], Use integrated 3D video system for observation without eyepieces (option), attach documentation / coobservation set-up.
 - ⇒ All necessary components are assembled on the surgical microscope and assistant's view.
 - ⇒ The tubes and eyepieces are correctly configured.
 - ⇒ Check that the surgical microscope and the tube are in a position convenient for you.
 - ⇒ The intraocular distance (pupillary distance) has been set.
 - ⇒ The eyecups have been set in such a way that you can see the full field of view.
 - ⇒ Check that the correct prescription has been set on the diopter scale.
 - ⇒ Check that image quality is correct throughout the entire magnification range.
 - ⇒ The objective lens and eyepieces are clean.
 - ⇒ The correct function of the other equipment (illumination system, video system, etc.) has been checked using the relevant Instructions for Uses.
 - ⇒ The surgical microscope can be positioned flawlessly with assembled accessories.
- ▶ Check handgrips on the surgical microscope.
 - ⇒ The buttons and rocker switches are equipped with the required functions.
 - ⇒ The buttons and rocker switches carry out the respective functions correctly.
- ▶ Check foot control panel and foot rocker switch.
 - ⇒ The buttons and rocker switches are equipped with the required functions.
 - ⇒ The buttons and rocker switches carry out the respective functions correctly.
- ▶ Check zoom [178].
 - ⇒ Start value and speed are set correctly.

- ▶ Check focus [178].
 - ⇒ Start value, speed and autofocus are set correctly.
- ▶ Check SmartLinks [179].
 - ⇒ The required SmartLinks are activated.
- ▶ Check displays.
 - ⇒ The inputs and outputs from/to external systems are correctly configured.
 - ⇒ The required images are displayed on the monitors/ external systems as required.
- ▶ Check image quality of the monitors.
 - ⇒ White balance was completed [196].
 - ⇒ Color saturation and hue are set correctly [195].
- ▶ Check image quality of the external systems, observe the manufacturer's specifications while doing so.
- ▶ If you wish to record [196] or stream [216] video data, check the respective settings (resolution, 2D/3D, parallel recording etc.).
- ▶ Check illumination [175].
 - ⇒ The illumination is set correctly.
 - ⇒ The remaining service hours are sufficient for the duration of the operation.
- ▶ Check user profile [201] and patient profile [205].
 - ⇒ The correct user is activated and correctly configured.
 - ⇒ The patient has been created correctly.

6 Operation

6.1 Safety during operation

WARNING!

Contamination with pathogenic germs!

If the device is cleaned using unsuitable methods or cleaning agents, pathogenic germs can accumulate on the device.

- ▶ Make sure that the device is operated only in a sterile manner and with sterile accessories.
- ▶ Make sure that the device is operated and cleaned only by instructed personnel.
- ▶ Use only cleaning agents which are suitable for cleaning sterile equipment when cleaning the device.

CAUTION!

The user interface of the touchscreen is unsterile!

The patient may be contaminated if sterile persons touch the touchscreen and then touch the sterile area.

- ▶ Ensure that the touchscreen is operated only by unsterile persons who have no contact with the sterile area.
- ▶ Ensure that sterile persons who have contact with the patient operate the touchscreen only via commercial, sterile medical applicators.

CAUTION!

Risk of injury to eyes caused by the light source.

If light enters the patient's eyes during surgery in the facial area, there is a risk of damaging the retina.

- ▶ Make sure that no light enters the patient's eyes during operations in the facial area.
- ▶ Do not use the device to perform ophthalmological interventions on the patient.

CAUTION!

Risk of injury to eyes caused by laser beams.

If laser beams enter the patients' eyes while the Focus or Autofocus are switched on, there is a danger of damage to the retina.

- ▶ During operations in the facial area, make sure that no laser beams from Focus or Autofocus enter the patient's eyes.
- ▶ While Focus or Autofocus are switched on, do not aim the surgical microscope at the patients' or any other persons' eyes.

⚠ CAUTION!

Risk of burns due to high light intensity.

The interaction of heat and antimicrobial substances in incision foils may lead to an increased reaction of the patient to these substances.

- ▶ Avoid using adhesive foil/incision foil to cover the patient wherever possible.
- ▶ Avoid high light intensity outside of the OR area.
- ▶ Avoid long treatment periods.

⚠ CAUTION!

Risk of burns due to high light intensity.

Using dry drapes to cover the patient increases the risk of 2nd degree burns due to an excessive light intensity.

- ▶ Avoid drying out of the drapes which are used to cover the patient.
- ▶ Make sure that the drapes always remain moist.
- ▶ Avoid high light intensity outside of the OR area.
- ▶ Avoid long treatment periods.

⚠ CAUTION!

Risk of burns due to high light intensity.

Using dry drapes to cover the patient increases the risk of 2nd degree burns due to an excessive light intensity when the device is switched on and left unattended.

- ▶ Never leave the device aimed at the patient while it is switched on.
- ▶ Never leave the device unattended for a longer period of time.

⚠ CAUTION!

Xenon lamp failure caused by ageing of the lamp.

If the xenon lamp is used beyond its maximum life of 500 service hours, it may fail suddenly and interrupt operation.

- ▶ Please note the remaining service hours displayed on the touch-screen.
- ▶ Once the current xenon lamp has reached 500 service hours, activate the backup xenon lamp.
- ▶ Replace the xenon lamp when its operating time of 500 service hours has expired.

⚠ CAUTION!

Damages due to unintentionally touching the touchscreen.

Touching the screen in this case may unintentionally open menus or change settings.

- ▶ Do not touch the touchscreen unintentionally if the system is switched on.
- ▶ Allow only trained personnel to use the system.

 CAUTION!

Risk of injury caused by electrical voltage.

Contact with plug connector contacts may cause electric shock.

- ▶ Never touch the contacts on plug connectors during contact with the patient.

 CAUTION!

Risks resulting from the drape glass cover being dirty

If the glass cover is dirty, this may interfere with the device's optical characteristics and prevent the device's functions (e.g. autofocus) from being used safely.

- ▶ Keep the drape glass cover free of dirt at all times.

6.2 Software configuration

6.2.1 Touchscreen

 CAUTION!

Do not touch the touchscreen if it has failed!

There is a possibility that only the illumination of the display has failed. In this case, you would access menus or change values inadvertently.

- ▶ Do not touch the surface of the touchscreen if the device is switched on and its screen is dark.

6.2.2 Configure microscope

Tap on  Settings →  Microscope.

Prior to use, check the illumination and make sure that the illumination intensity required for the application is available.

6.2.2.1 Configuring the illumination

You can assign the following illumination settings user specifically in the "Light" menu:

- Light on/off
- Additional light on/off
- Set light intensity


Figure 88: "Settings Microscope" menu, "Light" submenu

If a high magnification is used, the diameter of the field of view and the light intensity at the surgeon's eye both decrease whereas the light intensity in the surgical field remains the same. Some components can further reinforce this effect:

- Eyepieces with a high magnification
- 3-Step magnification changer
- Foldable tube f170/f260 with tube magnification (PROMAG function)

In digital operation, the video image can show more noise (autogain) or become darker (manual gain) at high magnification.

The "Focus Light Link" function limits the maximum adjustable light intensity in the focal plane depending on the working distance. Possible tissue damage due to an unintentional high light intensity can thus be avoided. This function should only be deactivated if a higher light intensity is required for the current application. Therefore, please note the recommendations for preventing burns.

Action

1. Tap on Settings → Microscope → Light.
2. Whenever the illumination is switched on, you also can switch on the additional light:
→ Additional light switch → activation .
3. **CAUTION! Tissue damage due to high light intensity!**
When working at maximum magnification, pay attention to the set light intensity to prevent burns, especially of the surrounding tissue.
Set the desired light intensity with the [Light Intensity] slide control.
The light intensity is influenced by the "Focus Light Link" function (if activated), so that the light intensity which you have set can be limited depending on the working distance.

4. Set the required warning threshold level on the [Light Warning Threshold] slide control.
 - ⇒ The value set is displayed in the "Light" menu.
If the light warning threshold is exceeded, this is displayed by an orange color in the live menu, in the main menu bar and on the Light Intensity slide control.
5. Select the speed of light intensity change with the buttons on the hand grip and the foot control panel, i.e. [High] or [Low].
6. In the "SmartLinks" menu you can activate or deactivate the "Focus Light Link" function.
 - ⇒ Activated: Influences the maximum adjustable light intensity and limits it depending on the working distance used.
7. In the "SmartLinks" menu you can activate or deactivate the "Auto Brightness" function.
 - ⇒ Activated: Automatic control of the light intensity for constant image brightness in the eyepiece.

* In China the "Focus Light Link" function always remains activated and cannot be deactivated due to national regulations.

6.2.2.2 Configure tube

Prerequisite

- ☒ The microscope is configured as a system for hybrid (optical-digital) visualization.

Action

- ☒ The system mode is set at "Optical".

1. Tap on  Settings →  Microscope → Tube [79].
2. Tap on the focal length of tube in use: 170 mm for tiltable tube and foldable tube or 260 mm for foldable tube with magnification.
3. Tap on the ocular magnification of the eyepieces used: 10x or 12.5x.
4. Tap on the magnification of the optional 3-position magnification changer.
5. Tap on the co-observation selection:
Opposite image outputs: Face-to-face for tiltable tube or foldable tube.
Lateral image outputs: Right and left for stereo co-observation module or photo adapter for external camera.

6.2.2.3 Configure focus

Action

1. Tap on  Settings →  Microscope → Focus.
2. Set the desired value with the [Start Value] slide control.
 - ⇒ The microscope starts with this working distance after it is switched on.
3. Set the desired value with the [Speed] slide control.
 - ⇒ Set speed of the motorized fine focus adjustment.
4. Set the desired value on the slide control [Autofocus speed] (option).
 - ⇒ Adjusting the Autofocus speed
5. Tap on the [Autofocus] button (option).
 - ⇒ On: Autofocus is automatically activated each time the brakes are closed.
 - ⇒ Off: Autofocus is triggered via a configurable button on the handgrip or the foot control panel.
The autofocus is not available during a few specific applications, e.g. when using a micromanipulator (Focus Lock function). A corresponding acoustic signal is then emitted.
6. Tap on the [laser focusing aid] button (option).
 - ⇒ Laser Focus Aid on: The focusing laser is switched on when the brakes are released or the motorized fine focus adjustment is actuated.
 - ⇒ Laser Focus Aid off: The focusing laser is switched off.
7. Tap on the [Focus stop] button. For the application of an optional micromanipulator with a fixed working distance.
 - ⇒ Focus Lock on: Autofocus inoperative. Focus rocker switch on hand grip/foot control panel out of operation.
 - ⇒ Focus Lock off: Autofocus functioning. Focus rocker switch on hand grip/foot control panel active.
8. Select the depth of field. Tap on the desired selection field in the "Focus Depth" selection field.
 - ⇒ High: high depth of field - less light - lower image resolution
 - ⇒ Low: low depth of field - more light - higher image resolution

The "Focus Depth" field is not displayed when the BLUE 400 and YELLOW 560 options are activated for the system.

6.2.2.4 Configuring zoom

Action

1. Tap on  Settings →  Microscope → Zoom.
2. Set the desired value on the [Start value] slide control.
3. Set the desired value on the [Speed] slide control.

6.2.2.5 Configuring SmartLinks

Action

1. Tap on  Settings →  Microscope.
⇒ In the "Microscope" menu you can activate or deactivate four SmartLinks.
2. Deactivate the "Focus Light Link"* function (this function is always activated by default).
⇒ A warning is displayed indicating that tissue damage may result if the "Focus Light Link" is deactivated.
3. If you really want to deactivate the "Focus Light Link" function, tap on the [Yes] button.
⇒ The "Focus Light Link" function no longer can limit the maximum adjustable light intensity in the focal plane (depending on the working distance). Possible tissue damage due to an unintentional high light intensity can thus occur.
This function should only be deactivated if a higher light intensity is required for the current application!
4. Therefore, please note the recommendations for preventing burns [18].
5. Activate the "Focus Light Link" function.
⇒ The maximum light intensity is limited for the selected working distance (zoom).
⇒ A high light intensity over a longer period of time can damage tissue and cause burns in the area of the illuminated surgical field. The "Focus Light Link" function reduces this risk and should always remain activated.
6. Activate the "Focus Zoom Link" function.
⇒ The focusing speed automatically adapts to the magnification. The preselected focusing speed is automatically reduced when a higher magnification is used.
7. Activate the "XY Zoom Link" function.
⇒ The motorized XY travel speed automatically adapts to the magnification. The preselected travel speed is automatically reduced when a higher magnification is used.
8. Activate the "Auto Brightness" function.
⇒ The light intensity automatically adapts to constant brightness in the eyepiece depending on the working distance and magnification.
9. Close the "Settings" menu by tapping on the  button.

* In China the "Focus Light Link" function always remains activated and cannot be deactivated due to national regulations.

6.2.3 Configuring the button assignment of the hand grips, foot control panel and rocker foot switch

6.2.3.1 Configuring the handgrips

In the "Handgrips" menu, you can assign user-specific functions to the buttons.

Action

1. Tap on Settings → Operation → Hand Grips.

2. To assign a function to a button: Tap on the corresponding letter e.g. [A].
⇒ A field with the available functions appears in the selection field on the right side.
3. Tap on the function to be assigned to the button.
⇒ The function is accepted and displayed next to the button symbol.
4. You can swap the zoom and focus functions of the two rocker switches and change the respective direction (+/-).
5. Tap on the rocker switch symbols [E-F], [G-H].
⇒ The available functions appear in the selection field on the right side.
6. Tap on the function to be assigned to the rocker switches.
⇒ The function is accepted and displayed next to the button symbol.

One of the configurable handgrip buttons also can be configured to change between the motorized movement modes ("Stand" XY Adjustment Mode menu).

7. To assign a function to a button: Tap on the corresponding letter e.g. [A].
8. Tap on the [XY Mode] function in the selection field, see XY mode movement modes [184].

6.2.3.2 Configuring the brake release buttons

The brake release buttons are located on the back of the handgrips

- The lower button [AB] is not configurable. When pressed, it releases all microscope and stand axes for positioning.
- The upper button [SB] (stand brakes) can be used to enable various movement modes, depending on the system configuration.

Action

1. Tap on  Settings →  Operation → Hand Grips.

2. To assign a function to the [SB] button: Tap on an [SB] button symbol.
 - ⇒ The available functions appear in the selection field on the right side.
3. In the selection field, tap on a function that you wish to assign to the button (e.g. SB: release stand brakes).
 - ⇒ The function is accepted and displayed next to the button symbol.

HINT: Test the button assignment and the functions of the brake release buttons before every application and without a patient.

6.2.3.3 Configuring the foot control panel (FCP)

In the "FCP" menu, you can assign user-specific functions to the buttons.

Action

1. Tap on  Settings →  Operation → FCP .
2. To assign a function to a button: Tap on a button symbol, e.g. [A].
 - ⇒ The available functions appear in the selection field on the right side.
3. Tap on a function in the selection field to be assigned to the button.
 - ⇒ The function is accepted and displayed next to the button symbol.
4. The zoom and focus functions can be swapped or the direction (+/-) can be changed.

5. Tap on the rocker switch symbols [I-G], [J-H].
⇒ The available functions appear in the selection field on the right side.
6. Tap on a function in the selection field to be assigned to the rocker switches.
⇒ The function is accepted and displayed next to the button symbol.

Furthermore, one of the configurable FCP buttons also can be used to change between motorized movement modes ("Stand" menu XY Mode).

7. To assign a function to a button: Tap on the corresponding letter e.g. [A].
8. In the selection field, tap on the function:[Mode XY adjustment].
9. Check the button assignment and the functions of the foot control panel before every application and without a patient.

6.2.3.4 Configure rocker foot switch

Action

1. Tap on  Settings →  Operation →  Rocker foot switch.
2. To assign a function to the rocker foot switch: Tap on button symbol [A] or [B]
⇒ The available functions appear in the selection field on the right side.
3. Tap on a function in the selection field to be assigned to the button.
⇒ The function is accepted and displayed next to the button symbol.

When the following functions are selected for a rocker foot switch, the second rocker foot switch is automatically assigned:


- Light: Light brighter (top button) / Light darker (bottom button)
 - Zoom: + (top button) / - (bottom button)
 - Focus: + (top button) / - (bottom button)
4. Check the button assignment and the functions of the rocker foot switch before every application and without a patient.

6.2.3.5 Protect Settings


Each user can protect his hand grip and foot switch configuration against third-party access using a password (PIN). After unlocking, the settings secured in this way can only be changed again using a valid PIN.

Action

1. Tap on Settings → Operation.
2. Tap on the [PIN input] button to activate the lock .

3. Enter a 4-digit PIN.

4. To save the PIN number: tap the [Done] button.
 - ⇒ The switch in the [PIN input] button lights up blue .
 - In the "Operation" menu, this is displayed for all controls at the upper right by the symbol.
 - The PIN lock is jointly **active** for all controls.
5. To **Unlock** tap on the [PIN input] button .
6. Enter your PIN.

7. Tap on the [Done] button.

- ⇒ The switch in the [PIN input] button lights up gray  .
The PIN lock is jointly **deactivated** for all controls.
In the "Operation" menu, this is shown in all controls, top right, by the symbol .

8. Close the "Settings" menu by tapping on the  button.

6.2.4 Configure stand

Action

1. Tap on  Settings →  Stand.

6.2.4.1 XY Mode movement modes

- Various brakes are activated in the manual movement modes. In the process, the microscope is guided via the hand grips by pushing the "AB" (all microscope and all stand brakes activated always) and the "SB" buttons.

The manual movement modes can be configured in  Settings →  Operation → Hand grips → SB button.

- In the motorized movement modes, the microscope movements are carried out using motor function by operating the joystick on the hand grip or FCP.

The motorized movement modes can be configured in  Settings →  Operation → Hand grip or FCP. The speeds of the motorized movement modes can be set via the "Stand" menu or in the main menu list on the touchscreen (see image).

One of the configurable handgrip or FCP buttons also can be configured to change between the motorized movement modes. ( Settings →  Operation → Hand grips → XY adjustment mode).

6.2.4.1.1 Manual movement modes

- Release "Stand Brakes"

In this mode, the brakes of the stand axes are released. The microscope remains fixed in its axes and can be moved freely about stand axes 1, 2 and 3 when the "SB" button is pressed on the handgrip.

- Release "microscope brakes"

In this mode, the brakes of the microscope are released with fixed stand axes. The microscope can be moved freely about axes 4, 5 and 6 when the "SB" button is pressed on the handgrip.


Figure 89: Overview of rotary axes: Stand axes: 1-3, microscope axes: 4-6

■ "PointLock"

If the "SB" button on the handgrip is pressed, the microscope can be moved and repositioned while the focal point remains in the center of the field of vision.

The working distance can be adjusted within the available focal length range of 200 ... 625 mm while the microscope is moved about the object to be observed. If you change the working distance while pivoting, the focus is automatically reset and the object in the center of the field of view remains sharply focused.

6.2.4.1.1.1 Configure manual movement modes

Action

1. Tap on ⚙ Settings → 🖱 Operation → Hand Grips.
2. Tap on the [SB] button. This is the green brake activation button at the top of the hand grip's reverse side.
 - ⇒ The assignable functions appear on the right-hand side:
 - "Stand brakes"
 - "Microscope brakes"
 - "PointLock"
3. Tap on the desired function.
 - ⇒ The assigned functions appear alongside the [SB] button.

4. Close the "Settings" menu by tapping on the  button.
5. On the hand grip, press the brake activation [SB] button with the required function and continue to hold it pressed.
 - ⇒ "Stand brakes": Only the stand axes 1, 2, and 3 can now be freely moved manually.
 - ⇒ "Microscope brakes": Only the microscope axes 4, 5, and 6 can now be freely moved manually.
 - ⇒ "PointLock": Manual movement of the microscope. In this movement mode the working distance can also be changed, the middle of the visual field remains automatically in focus and well defined (in the area of the working distance from 200...625 mm).

6.2.4.1.2 Motorized movement modes

6.2.4.1.2.1 Microscope movement mode

The microscope can be moved by motor in the three axes of its suspension, 4, 5 and 6. Precise motorized XY movement of the focal point is thus enabled by tilting and swiveling the microscope with a fixed working distance. The stand does not move in the process.

Movement can be triggered via the joystick on the handgrip or the FCP if the corresponding XY movement mode was selected in the "Stand" menu. The movement speed is individually adjustable.

Prerequisite

- ☒ The microscope is positioned and focused above the desired surgical field


Action

1. Tap on  Settings →  Stand.
2. Tap on the [Microscope] button in the "XY Mode" field.
3. Set the desired speed with the [Microscope] slide control in the "Speed settings (XY)" field.
 - ⇒ The assigned speed is displayed on the respective slide control.
4. Close the "Settings" menu by tapping on the  button.

6.2.4.1.2.2 Stand movement mode

The microscope can be moved by motor in the three axes 4, 5 and 6 in the XY focal plane without tilting or swiveling; the horizontal alignment of the eyepieces always remains unchanged. This enables precise motorized XY movement of the focal point at a fixed working distance.

Movement can be triggered via the joystick on the handgrip or the FCP if the corresponding XY movement mode was selected in the "Stand" menu. The movement speed is individually adjustable.

Prerequisite

- ☒ The microscope is positioned and focused above the desired surgical field


Action

- ▶ Tap on ⚙ Settings → 🏠 Stand.
- ▶ Tap on the [Stand] button in the "XY Mode" field.
- ▶ Set the desired speed with the [Stand] slide control in the "Speed settings (XY)" field.
 - ⇒ The assigned speed is displayed on the respective slide control.
- ▶ Close the "Settings" menu by tapping on the 🏠 button.

6.2.4.1.2.3 PointLock movement mode

This motorized movement function enables time-saving and precise movement of the microscope, during which the microscope remains focused on the focal point (center of the field of view). Movement can be triggered via the joystick on the handgrip or the FCP if the corresponding XY movement mode was selected in the "Stand" menu. The movement speed is individually adjustable. In the Motorized PointLock mode, you can use the joystick on the

handgrip or the FCP to motorically pivot the microscope in the same direction as the joystick (CCW/CW/forward/back – in 90° increments only!) and only with a fixed working distance!).

Prerequisite

- ☒ The microscope is positioned and focused above the desired surgical field

Action

1. Tap on Settings → Stand.
2. Tap on the [PointLock] button in the "XY Mode" field.
3. Set the desired speed with the [PointLock] slide control in the "Speed settings (XY)" field.
You also can set this speed via the main menu bar on the touchscreen.
⇒ The assigned speed is displayed on the respective slide control.
4. Close the "Settings" menu by tapping on the button.

6.2.4.1.2.4 Configuring the handgrip "XY Mode"

One of the configurable handgrip or FCP buttons also can be configured to change between the motorized movement modes.

Action

1. Tap on Settings → Operation → Handgrips or FCP.
2. To assign a function to a button of a handgrip or the FCP: Tap on the corresponding letter of the button e.g. [A].
⇒ In the "Handgrips" or "FCP" menu, a selection field with the available functions appears on the right side.
3. Tap on [XY Mode] in the selection field.
⇒ The function is accepted and displayed next to the button symbol.

Check the button assignment and the functions of the handgrips or the foot control panel before every application and without a patient.

6.2.4.1.2.5 Parking position

CAUTION!

Risk of crushing!

Fingers may be crushed between the vertical arm and the horizontal arm.

- ▶ The parking position may only be used after the application (without patient) or for system transportation!
- ▶ Never touch the area between the vertical arm and the horizontal arm while they are moving to the park position.

Action

1. Bring the monitor (optionally: both monitors) into the transit position.

2. Tap on  Settings →  Stand → Move to park position → Start.
3. To activate the travel, press the XY joystick button on the right handgrip to the right  and hold it down until the park position has been reached. The joystick on the FCP also can be used to activate the travel.
4. The device now automatically travels to the park position.

- ⇒ When the park position has been reached, this is acknowledged by a signal tone.

6.2.4.1.2.6 Drape position

The drape position is saved as a factory setting, however, can be changed arbitrarily.

CAUTION!

Risk of crushing!

Fingers may be crushed between the vertical arm and the horizontal arm.

- ▶ The drape position may be used only without a patient!
- ▶ Never touch the area between the vertical arm and the horizontal arm while you are moving to the drape position.

Action

1. Tap on  Settings →  Stand.
2. Tap on the [Start] button in the "Start Drape Position" field.
3. To activate the travel, press the XY joystick button on the right handgrip to the right  and hold it down until the drape position has been reached.
The joystick on the FCP also can be used to activate the travel.
 - ⇒ The device then automatically travels to the drape position.
 - ⇒ When a beep sounds, the drape position has been reached.

Determining a new drape position

If the factory set drape position is unfavorable for you, you can define and save your own drape position.

4. Move the device to a position that is favorable and comfortable for you.
5. Tap on  Settings →  Stand.
6. Tap on the [OK] button in the "Save Current Position as Drape Position" field.
 - ⇒ The current position of the device is saved as the new drape position.
7. Now you can move the device to the drape position as usual.

6.2.4.2 Position memory

CAUTION!

Tissue damage due to incorrect movement of the surgical microscope!

If the system or patient is moved or displaced in the meantime, it is longer possible to move to the saved system position safely.

- ▶ Do not move the system or the patient if you wish to use position saving.

You can intraoperatively save the current position of the microscope, its orientation to the object, the working distance and the magnification at any time during a surgical procedure via a preconfigured handgrip or FCP button or recall these parameters again directly via the monitor.

You no longer can move to the saved positions in the following cases:

- The device was restarted (the reference coordinate system of the stand was reinitialized)
- The power supply was interrupted (the reference coordinate system of the stand was reinitialized)
- The parking brakes of the device were released and the device was moved
- The patient was moved

6.2.4.2.1 Configuring the position memory


Action

1. Tap on Settings →  Operation → Handgrips or FCP.
2. In order to assign a hand grip or FCP button to position saving:
Tap on the corresponding letters of the button e.g. [A].
⇒ In the "Handgrips" or "FCP" menu, a selection field with the available functions appears on the right side.
3. In the selection field, tap on [Position Memory].
⇒ The function is accepted and displayed next to the button symbol.

6.2.4.2.2 Save Position

Action

1. Tap on the hand grip or FCP button configured with "Position Memory".
⇒ A miniature image of the photograph taken is briefly displayed on the screen [Photos] on the touchscreen. The photo is saved with the current microscope position, working distance and magnification in the patient data with the "Position" symbol .
2. You can also use the  button on the touchscreen in the main menu list below.
⇒ The "Position Memory" menu is shown on the bottom right-hand side of the touchscreen.

3. Tap on [Save position].

6.2.4.2.3 Start Position

Action

1. Press and hold the position button for at least 2 seconds. A selection menu will open on the touchscreen showing the positions currently saved.
2. Press and release the joystick button on the right hand grip or on the FCP to navigate to the desired saved position.
3. The active position is highlighted in blue.
4. Press the joystick button on the hand grip or the FCP to the right and hold it down until the position has been reached.
 - ⇒ When the position has been reached, this is acknowledged by a signal tone.

6.2.4.2.4 Set positioning speed

Action

1. Tap Settings → Stand → scroll down in the stand menu.
2. Set the desired positioning speed on the slide control.
3. Close the "Settings" menu by tapping on the button.

6.2.5 Configure displays

Display on external monitors, in INFRARED 800 option on the touchscreen

Action

1. Tap on Settings → Displays.
2. In this menu, set which display should be shown on external monitors at what level of quality.
3. Close the "Settings" menu by tapping on the button.

6.2.6 Configure audio

Action

1. Tap on Settings → Audio.
2. Internal microphone: → Switch on "Use switch for → video recording" .
3. Internal Microphone: → Switch on "Forward switch to line-out" → .
4. External microphone: → Switch on "Use switch for video recording" → .
5. External microphone: → Switch on "Forward switch to → line-out" .

6. External microphone: → Switch on "Forward switch to internal loudspeaker" →  →. Set the volume with the controller.
7. Volume: Video Playback: → Set the volume on the slide control.
8. Volume: System Sounds: → Set the volume on the slide control.
9. Close the "Settings" menu by tapping on the  button.

6.2.7 Configuring Photo

Action

1. Tap on  Settings →  Photo.
2. Select whether the photos are to be taken with the internal camera or with the external DSLR camera (option) in the "Recording" field.
3. Select the desired format for saving your photos on the device in the "Storage format (internal camera)" field.
4. If you would like to take photos with an external camera, activate the slide switch in the "Light intensity for external camera" field.
 - ⇒ The switch illuminates blue .
 - ⇒ The light intensity of the xenon illumination is briefly set to 100% when the external DSLR camera is released.
5. Close the "Settings" menu by tapping on the  button.

6.2.8 Configuring the video camera

All settings in this menu are displayed in the preview image.

Action

1. Tap on  Settings →  Video Camera.
2. Set the parameters of the video camera as desired.

6.2.8.1 Setting the brightness

You can set the exposure for the video recording via the shutter control and via the direct exposure control:

- Automatic shutter control: The video camera automatically controls the exposure time and keeps the set brightness constant
- Manual shutter control: You can set the exposure time manually
- "Spot" exposure control: The exposure is measured in a very small area in the image center. This setting is suitable for working in narrow channels.

- "Full" exposure control: The exposure is measured across the entire video image. This setting is suitable for fully illuminated surgical fields or strong local reflections.
- Auto Brightness: You can set the desired exposure time manually; the video camera automatically readjusts the exposure time

Action

1. Select the "Auto" or "Manual" function in the "Shutter Control" field.
2. Select the function "Spot" or "Full" in the "Exposure Control" field.
3. Set the required brightness with the slide control in the "Auto Brightness" field.

6.2.8.2 Setting the color saturation and hue

You can set the hue and color saturation (chroma) of the integrated camera.

Action

1. Set the desired value with the slide control in the "Color Saturation" field.
 - ⇒ Higher values produce more intensive colors.
2. Set the desired value with the "Hue" slide control.
 - ⇒ A hue of ± 0 usually provides the most natural color impression.

6.2.8.3 Setting the black level offset

The black level offset indicates how bright the darkest part of the image is.

Action

- Set the desired value with the slide control in the "Black Level Offset" field.
 - ⇒ The smaller this value is, the darker the black will be.

6.2.8.4 Setting the red level, blue level and peak/average

You can set the coloring of the integrated camera via the red level and the blue level. You can influence the characteristic of the automatic exposure via the peak/average level.

Action

- Set the desired values in the "Red Fine Adjustment", "Blue Fine Adjustment" and "Exposure Control Peak/Average" fields using the respective slide control in each case.
 - ⇒ Values of ± 0 usually provide the most natural color impression.

6.2.8.5 White Balance

You can perform a white balance for the integrated camera. Here the device adjusts the signal of the integrated camera in such a way that white areas in the surgical field also appear white on the monitor.

The white balance must not be performed during a fluorescence application in order to ensure that the camera parameters are not falsified. Furthermore, the influence of extraneous light should be minimized or avoided altogether (e.g. direct sunlight passing through a window or direct illumination from a surgical ceiling lamp).

Prerequisite

- ☒ The illumination is switched on on the device
- ☒ The device is not in the fluorescence mode
- ☒ The surgical field is shielded from extraneous light

Action

1. Aim surgical microscope at a white object.
2. Focus on the white object.
 - ⇒ The field of view is sharply discernible.
3. Start the white balance. Tap on the [White Balance] button.
 - ⇒ The following prompt appears: "Please place a sheet of white paper underneath the microscope, focus on it and press "Next!"
4. To start the white balance: Tap on the [Next] button. To cancel the process: Tap on the [Cancel] button.
 - ⇒ The message "Please wait - white balance in progress!" appears. After a successful white balance, the message "White balance completed successfully" appears. If it was not successful, the message "White balance failed" is displayed.
5. Following a successful white balance, tap on the [Close] button.
 - ⇒ The white balance settings are saved user-dependently.
6. Close the "Settings" menu by tapping on the  button.

6.2.9 Configuring the video recording

Action

1. Tap on  Settings  → Video Recording.
2. If the recorded video data are to be saved only on the internal hard drive, tap on "None" in the "Parallel Recording" field.
3. If the recorded video data also are to be (parallel) saved to a connected USB medium, tap on "USB" in the "Parallel Recording" field.
4. If the recorded video data also are to be (parallel) saved to a network folder, tap on "Network" in the "Parallel Recording" field.

5. Switch on "3D Recording" (option) →  if the video data are to be recorded in 3D format.
6. Switch on "SmartRecording" →  if you would like to record HD video data extending into the past for a period of up to five minutes during an operation.
 - ⇒ If "SmartRecording" is activated, HD video data of the operation are recorded in a buffer. The oldest video data in the buffer are cleared in such a way that the last five minutes always remain stored in the buffer. You can start and stop the recording the HD video data from the buffer at any time.
7. Switch on "Video Streaming" →  if you would like to live stream video recordings to a streaming client via LAN, WLAN or hotspot.
 - ⇒ The web address of the device which must be entered as the Streaming Client is displayed in the "Web Address for Video Streaming" field.
8. If you would like to transmit your video data for video streaming in high definition, tap on "HD Video" in the "Video Streaming" field.
 - ⇒ The video files are transmitted in high definition.
9. If you would like to transmit your video data for video streaming in low resolution, tap on "LowRes Quality" in the "Video Streaming" field.
 - ⇒ The video files will then be transmitted in low resolution. The quality of the transmitted video data is therefore lower than for recordings performed in HD quality.
 - ⇒ The data volume is much lower.
10. Close the "Settings" menu by tapping on the  button.

6.2.10 Configure fluorescence

All of the settings you select in the "Fluorescence" menu apply only to the currently active user.

Action

1. Tap on  Settings  → Fluorescence.
2. Select the desired fluorescence application in the "Fluorescence Mode" field.
 - Only activated fluorescence options are displayed.
3. BLUE 400: "Light Intensity (start value)" → Set the start value with the slide control.
4. YELLOW 560: "Light Intensity (start value)" → Set the start value with the slide control.

5. IR 800: Switch on "AutoGain" → .
⇒ The camera gain is controlled automatically during the recording phase.
6. IR 800: Switch off "AutoGain" → .
⇒ You can control the camera gain manually during the recording phase.
7. IR 800: Switch on "AutoDetection" (automatic fluorescence detection) → .
⇒ The "Replay" playback phase initially starts only when the fluorescence agent begins to flow into the tissue.
8. IR 800: Switch off "AutoDetection" (automatic fluorescence detection) → .
⇒ The "Replay" playback phase already starts at the beginning of the video recording, i.e. the black leader (where no fluorescence signal is present) also is played back.
9. IR 800: Switch on "Skip Setup" (preparation phase) → .
⇒ INFRARED 800 starts immediately in the Recording Mode. You cannot perform any focus or zoom settings.
10. IR 800: Switch off "Skip Setup" (preparation phase) → .
⇒ You can set the focus and zoom parameters and select the image section; then the recording phase starts.
11. IR 800: Switch on "AutoZoom" (default setting: off) → .
⇒ The system starts the recording phase with the magnification set by you in the "AutoZoom (start value)" field.
12. IR 800: Switch off "AutoZoom" → .
⇒ You can set the total magnification outside of the preconfigured range.
13. IR 800: "AutoZoom (start value)" → Set the desired start value with the slide control.
14. IR 800: "Gain (Start value)" → Set the desired start value on the slide control.
15. IR 800: "Short Replay [Number]" → Set the desired repetition rate for the short replays with the slide control. You can set the repetition rate to a value ranging from 0 to 10.
⇒ Short replays are played back successively until the repetition rate set by you is reached or the Fluorescence button on the handgrip or foot control panel is pressed.
16. IR 800: "Short Replay [Duration]" → Set the duration for the short replays with the slide control.

17. IR 800: "Long Replay [Number]" → Set the desired repetition rate for the long replays with the slide control.
You can set the repetition rate to a value ranging from 0 to 10.
⇒ Long replays are played back successively until the repetition rate set by you is reached or the Fluorescence button on the handgrip or foot control panel is pressed.
18. IR 800: "MultiVision overlay" → set whether the infrared image should also be shown during recording as MultiVision display () or shown on the monitor only ().
19. IR 800: "Display on the monitor" → set whether the infrared image should be shown on the touchscreen during recording.
20. FLOW 800: "MultiVision Data Injection" → set whether the FLOW 800 evaluation is also shown as a MultiVision Data Injection () or shown on the monitor only ().
21. Close the "Settings" menu by tapping on the  button.

6.2.11 Configure MultiVision

In this menu, you can set the parameters for data injection.

Action

1. Tap on  Settings →  Operation → Handgrips or FCP.
2. To assign the MultiVision function to a button of a handgrip or the FCP: Tap on the corresponding letter of the button e.g. [A].
⇒ In the "Handgrips" or "FCP" menu, a selection field with the available functions appears on the right side.
3. Tap on MultiVision in the selection field.
⇒ The function is accepted and displayed next to the button symbol.
4. Tap on  Settings ● → MultiVision → MultiVision switch → Switch-on  .
5. Select the MultiVision source. To do this, tap on the desired selection field.
6. Select the type of status displays on the display.
7. Use the [Brightness] slide control to set the desired brightness on the display.

8. If you would like to show MultiVision demo videos, tap on → Show MultiVision Demo Videos → Switch on .

- ⇒ The menu with MultiVision demo videos is displayed.
- ⇒ The "Overlay Mode (no Shutter)" field is displayed and the overlay mode is automatically switched on.
You can display the MultiVision demo videos with the shutter open or closed (with or without the overlay mode).
9. If you would like to view MultiVision demo videos with the shutter closed, tap on → Overlay Mode (no Shutter) → Switch off .
10. Tap on one of the video thumbnails.
⇒ The selected MultiVision demo video is shown.
11. Before you leave the MultiVision demo mode, make sure that the shutter is open or open the shutter: Tap on → Overlay Mode (no Shutter) → Switch on .
12. To activate the automatic switch-off function for QEVO (option), tap on → QEVO → Automatic Switch-off .
13. To switch on the marker for QEVO (option), tap on → QEVO → Marker .
14. Select the video source. To do this, tap on the desired selection field.
15. Depending on the navigation system connected, tap on → "Navigation (standard)"  or "Navigation (extended)" .
16. Close the "Settings" menu by tapping on the  button.

6.3 Manage User

In the "User" menu, you can manage a maximum of 40 users. Each user may be assigned a user language and various user-specific settings.

- You can change the names of users. You cannot change the name of the default user.
- You can delete users. You cannot delete the active user or the default user.
- If the desired user is not displayed, scroll down on the touch-screen display.

6.3.1 Adding a user

When adding a new user, the default settings are firstly adopted, user-specific settings (labelled on the touchscreen with a green dot) [64] may be carried out by the user at any time.

Action

1. Tap on the [User] button in the status bar.
2. Tap on the button .
3. Enter: name, language, interpupillary distance, diopters of the left/right eye, marking rings of the adjustable eyecup for the new user.
4. Tap on the [OK] button .
5. In order to select the new user as an active user, tap on the desired user name in the "Select user" menu.
 - ⇒ The current user is highlighted in blue.
 - ⇒ The new user is now highlighted in blue and is shown as an active user in the status list.
6. Close the menu by tapping on the button .

6.3.2 Configuration of users

In the status bar of the "Configure User" menu , you can perform the following user profile settings for the current user:

- Editing the user name
- Language selection
- Add comment
- Enter the interpupillary distance
- Enter left eye diopters
- Enter right eye diopters
- Marking ring on the adjustable eyecup

If the entered values (in the case of interpupillary distance, diopters of left and right eye) are greater than adjustable values, a note on the permissible values appears at the top of the touchscreen.

1. Tap on the [Configure User] button  in the status bar.
2. Tap on the respective field: [user name, language, comment ...].
3. Enter the information for the user profile.
4. Tap on the [Save] button.
 - ⇒ The user profile settings are now saved for the current user.

6.3.3 Export user settings to a USB storage medium.

Every user (incl. default users) can export his own complete user settings to a USB storage device in order to transfer them to another KINEVO 900. Only the IT admin user can at once select and export all user data records to a USB storage device.

1. Select the user to be exported as the active user.
2. Tap on the [Configure User] button  in the status bar.
3. Connect the USB storage device to the USB port.
4. Tap on the [Export] button  to export the complete user settings to a USB storage device.

Action

Action

6.3.4 Importing the user settings from a USB storage device (e.g. from another KINEVO 900)

A logged in user (incl. Default_User) may import his user settings from a user data record from a USB storage medium. Only the IT Admin can import complete user data records, with these then being created as new users.

When a new user is created, it is always checked whether this name already exists. If so, an addition is applied to the newly added user name (e.g. "Ellen Edwards_2").

Action

1. Tap on the [Configure User] button  in the status bar.
2. Connect the USB storage device to the USB port.
3. Tap on the [Import] button  to import the complete user settings from the USB storage device.
 - ⇒ The list of user records stored on the USB storage device is displayed.
4. Select the user record to be imported.
 - ⇒ The selected user record is imported.
The selected user is displayed as the active user in the status bar.

6.3.5 Changing user

The IT admin password [149] must be entered after user selection in order to activate users with IT admin system rights.

Prerequisite

- ☒ No video recordings are active.

Action

1. Tap on the [User] button in the status bar.
 - ⇒ The "Select User" menu is displayed. The current user is highlighted in blue.
2. Tap on the user to be selected.
 - ⇒ The new user is displayed in the status bar. The "Change of user" menu is hidden.

6.3.6 Copy a user

Action

1. Tap on the [User] button in the status bar.
 - ⇒ The "Select User" menu is displayed.
All stored users are displayed. The user list is sorted alphabetically.
2. If a user is not displayed, scroll down the user list.
3. Tap on the [Configure User] button .
4. Mark the box in front of the user name.
 - ⇒ The desired user is highlighted in blue .
5. Tap on the [Copy] button  at the bottom left to copy the user record.
 - ⇒ The copy (with the suffix "_1" appended to the name) is displayed below the copied user.

6.3.7 Delete User

Action

1. Tap on the [User] button in the status bar.
 - ⇒ The "Select User" menu is displayed.
2. Tap on the [Configure User] button .
3. If the new user is not displayed, scroll down the user list.
4. Mark the box in front of the user name to be deleted.
To delete multiple users simultaneously, mark the boxes in front of the names of all users to be deleted.
 - ⇒ Your selection is confirmed by a  symbol.
5. Tap on the [Delete] button  to delete the user name(s) you have selected.
 - ⇒ The "Delete User" menu is displayed.
6. To delete the user(s): Tap on the [Yes] button.
 - ⇒ The user(s) has/have been deleted.
7. Tap on the live image on the touchscreen to close the menu.

6.4 Patient administration

You can manage patient data in the "Patient" menu. In the patient menu, you can save, export or delete images, videos or data.

6.4.1 Adding a patient

Action

1. Tap on the [Patient] button in the status bar.
⇒ The "Patient" menu is displayed.
2. Tap on the [Add Patient] button .
⇒ The "Create New Patient" menu is displayed.
3. Tap on the respective field: [First name, Last name, Gender, Date of birth, Patient ID, Comments ...].
⇒ The virtual keyboard is displayed.
4. Enter the desired patient data.
If the date is entered in the wrong format, a message will appear at the top of the touchscreen: "Wrong format. Use MM/DD/YYYY."
5. Tap on the [Done] button.
⇒ The patient data were saved for the current patient.
6. To close the "Patient" menu: Tap on the [Patient] button or on the  button in the status bar.

6.4.2 Changing patients

It is not possible to change patients during a video recording. Always terminate any running video recordings before changing patients.

Action

1. Tap on the [Patient] button in the status bar.
⇒ The "Patient" menu is displayed.
2. Tap on the new patient name.
⇒ The patient folder is displayed and created as the current patient.
3. To close the "Patient" menu: Tap on the [Patient] button or on the  button in the status bar.

6.4.3 Edit patient data

In the status bar of the "Edit patient data" menu , you can perform the following settings for the current patient:

- Change first name
- Change last name
- Add comment
- Gender

Action

- Date of birth
 - Patient ID
1. Tap on the [Edit patient data]  button in the status bar.
⇒ The "Edit patient data" menu is displayed.
 2. Tap on the respective field: [First name, Last name, Gender, Date of birth, Patient ID, Comments ...].
⇒ The virtual keyboard is displayed.
 3. Change the required patient information.
 4. Tap on the [Done] button.
⇒ The patient data were saved for the current patient.

6.4.4 Searching for a patient

Action

1. Tap on the [Patient] button in the status bar.
⇒ The "Patient" menu is displayed.
2. Tap on the [Search] button .
⇒ A search window is displayed.
3. Enter a search term: Patient name, starting letters, etc.
⇒ The patient or patient group being searched for is displayed.

6.4.5 Sorting patients

Action

1. Tap on the [Patient] button in the status bar.
⇒ The "Patient" menu is displayed.
2. Tap on the [Filter] button .
⇒ The sorting criteria [Date], [Name], [Size] are displayed.
3. Tap on the desired [Button].
⇒ The patient list is displayed sorted.
4. Tap on the [Sort] button  to display the sorted patient list in ascending or descending order (alphabetical order).
5. To close the "Patient" menu: Tap on the [Patient] button or on the  button in the status bar.

6.4.6 Export Patient Data

Note:


Patient data that are exported from a system with software version 1.6 or higher cannot be imported into systems that have a lower software version than 1.6.

Prerequisite

- ☑ When exporting to a USB storage medium, the target storage medium has enough free space; the USB symbol in the status bar must not be colored orange.

Action

1. Tap on the [Patient] button in the status bar.
⇒ The "Patient" menu is displayed.
2. Tap on the name of the patient whose data is to be saved. Mark the box in front of the patient name.
⇒ Your selection is confirmed by a  symbol.
3. Tap on the [Export Selected] button .
⇒ The "Export Patient Data" menu is displayed.

4. Select the desired export path: "USB", "Network" or "DICOM" button.

Note on export to a USB storage medium:


If your USB storage medium does not have sufficient storage space and you start an export, an "insufficient storage space" error message will initially appear during, and not at the beginning of, the export process.

5. Check the available space on your USB storage medium before starting an export.
6. Select the memory option:
DICOM format (for USB and network storage only)
 On: Images are stored on a USB medium in DICOM format.
 Off: Images are stored in native format (PNG, JPEG) in the predefined patient data structure.
Anonymized
 On: No user or patient names are stored in the EXIF image data.
 Off: The EXIF image data include surgeon and user names.
7. Tap on the [Start] button to export the patient data.
8. To close the Patient menu: Tap on the "Patient" button or on the  button in the status bar.
9. A successful export is acknowledged by a short status message.

6.4.7 Importing patient data from the USB port or network

Action

1. Tap on the [Patient] button in the status bar.
⇒ The "Patient" menu is displayed.
2. Tap on the [Import]  button.
3. Select the medium from which the patients are to be imported:
[USB] or [Network] button.


- ⇒ The "Import Patient Data" menu is displayed.
4. Mark the box in front of the name of the patient ☒ whose patient data you want to import.
Mark the [Select All] box if you want to import all of the displayed patients.
 5. Tap on the [Start] button to import the patient data.
⇒ A new patient record is always imported, i.e. the data are not added to the existing patient record.
⇒ If a new patient record is imported which has the same name as the existing patient record, there will be two patient records with the same name.
If necessary, change the new patient name (by adding the suffix "_1").
 6. To close the "Patient" menu: Tap on the [Patient] button or on the  button in the status bar.

6.4.8 Importing patient data from the DICOM Worklist

Action

1. Tap on the [Patient] button in the status bar.
⇒ The "Patient" menu is displayed.
2. Tap on the [Import]  button.
3. Tap on the [DICOM] button.
⇒ "Importing patient data from the DICOM Worklist" input mask is displayed.

Scheduled Date	Last Name	First Name	Patient ID	
18.07.2018				☐
19.07.2018				☐

4. If you know the patient's required data and would like to import it, enter this data in the input mask.
5. If you do not know the patient's required data, tap on the [Show Worklist] button.
⇒ All patients of the DICOM Worklist are displayed.
6. If you would like to view further details of individual patients, tap on the corresponding arrow key on the right side of the screen.
⇒ A detailed list of the respective patient data then opens.
7. Push this list to the right to return to the DICOM Worklist.
8. Mark the box in front of the name of the patient ☒ whose patient data you want to import.
Mark the box in the header of the DICOM Worklist if you want to import all of the displayed patients.

9. Tap on the [Import selection] button.
 - ⇒ The "Patient" menu is displayed.
 - ⇒ The imported patient(s) are assigned a "D" in the patient list.

 - ⇒ The patient name will now have an information symbol next to it in the status bar.

 - Tapping this symbol opens a detailed list of patient data.
10. To close the Patient menu: Tap on the "Patient" button or on the  button in the status bar.

6.4.9 Deleting patients

Action

1. Tap on the [Patient] button in the status bar.
 - ⇒ The "Patient" menu is displayed.
2. Tap on the patient to be deleted.

To delete multiple patients simultaneously, mark the boxes in front of the names of all patients to be deleted.

 - ⇒ Selection is confirmed by a  symbol.
3. Tap on the [Delete]  button to delete the selected patient.
 - ⇒ The "Delete Patient Data" menu is displayed.
4. To delete the patient(s):

First place a check mark  next to "Remove patient". Then tap on the [Yes] button.

 - ⇒ The patient or patients is/are then deleted.
5. To close the "Patient" menu: Tap on the [Patient] button or on the  button in the status bar.

6.5 Recording video files (option)

With the "Video Recording" option, you can record the live image of the integrated video camera during an operation. The recorded video file is respectively assigned to the current patient in the process. A video file can be recorded for up to 24 hours.

You can record video files in two different resolutions:

- HD quality
 - The video files are recorded in high definition
 - The data volume increases in proportion to the length of the video recording
 - Recording video data in HD quality is suitable for short operations

- Video recordings in HD quality which extend up to 5 minutes into the past can be performed by activating the "SmartRecording" function [214]
- LowRes quality (SD)
 - The video files will then be recorded in low resolution. The quality of the recorded video data is therefore lower than for recordings performed in HD quality.
 - The data volume is much lower.
 - Recording video data in low resolution quality is suitable for long operations.

Prerequisite

- ☑ The parameters in the "Video Camera" menu have been set
- ☑ The parameters of the "Video Recording" menu have been set [196]
- ☑ The "Video Recording" function has been assigned to a button on the handgrip or the foot control panel
- ☑ The internal hard drive has enough free space; the hard drive symbol in the status bar must not be colored orange

Action

1. Start the recording of the video files in high definition (HD) by tapping on the [Video Recording] button  on the screen or by actuating the preconfigured [Video Recording] button on the handgrip or on the foot control panel.
 - ⇒ On the screen, the [Video Recording] button changes from  to  and the time display indicates the runtime of the video recording.
 - ⇒ A blue bar with information on the video recording appears at the bottom of the screen:
On the left side: HD recording; on the right side: storage location, e.g. internal hard drive with specification of storage capacity in percentage.

 - ⇒ The blue bar with information on the video recording turns orange as soon as the free storage capacity of the internal hard drive is only 20%.
2. As soon as the blue bar with information on the video recording turns orange, stop the recording of the video files and make sure that there is sufficient storage space on the internal hard drive.
3. Stop the recording of the video files in high definition (HD) by tapping on the [Video Recording] button  on the screen or by actuating the preconfigured [Video Recording] button on the handgrip or on the foot control panel.

4. To start recording the video files in low resolution (LowRes/SD), open the additional menu "Extended operation" and tap on the [LowRes Recording] button .

⇒ On the screen, the [LowRes Recording] button changes

from  to  and the time display indicates the runtime of the video recording.

⇒ A blue bar with information on the video recording appears at the bottom of the screen:

On the left side: LowRes recording; on the right side: storage location, e.g. internal hard drive with specification of storage capacity in percentage.

⇒ The blue bar with information on the video recording turns orange as soon as the free storage capacity of the internal hard drive is only 20%.

5. As soon as the blue bar with information on the video recording turns orange, stop the recording of the video files and make sure that there is sufficient storage space on the internal hard drive.
6. Stop the recording of the video files in low resolution (LowRes/SD) by tapping on the [LowRes Recording] button  on the screen.

Hint

You can take additional photos or capture images of interesting points while recording HD video files. These photos are assigned a marker in the HD video files (photos are marked in video files only for HD videos and not for LowRes recordings). During the playback of recorded HD video files, small gray dots are displayed above the "Video playback" slide control at the locations where you have taken a photo.

7. Tap on the [Photo] button  on the screen or actuate the preconfigured [Photo] button on the handgrip or foot control panel.
- ⇒ An acoustic signal sounds.
- ⇒ The photo is then saved in the "Patient" menu under the current patient.
- ⇒ This location is furnished with a marker in the video file and marked by a small gray dot located above the playback bar.
8. Open the recorded video file.
9. Move the "Video playback" slide control to a location which is marked by a small gray dot.

6.6 Additional recording of video files on a USB storage medium or network drive (parallel recording)

⇒ You now can view the photo.

6.6 Additional recording of video files on a USB storage medium or network drive (parallel recording)

Parallel to recording a video on the internal hard drive of the device, you can record the video files on a connected USB storage medium or network drive.

6.6.1 Parallel recording of video files on a USB storage medium

Note: It is not possible to use USB sticks for parallel recording.

Prerequisite

- ☑ The USB storage medium uses at least USB 3.0 (e.g. the USB 3.0 hard drive supplied with the device).

Action

1. Connect a USB storage medium to the device.
⇒ The USB symbol and the amount of storage space occupied on the USB storage medium appear in the status bar on the upper right edge of the screen.
2. Make sure that enough storage space is available; i.e. the USB symbol must not be colored orange.
3. Tap on ⚙ Settings → Video Recording.
4. Activate the "USB" function in the "Parallel Recording" field.
5. Close the "Settings" menu by tapping on the ⬆ button.
6. Start the recording of the video files [210].

6.6.2 Parallel recording of video files on a network drive

Prerequisite

- ☑ The connection to the network drive (Network Storage) has been established and is stable [152]

Action

1. Tap on ⚙ Settings → Video Recording.
2. Activate the "Network" function in the "Parallel Recording" field.
3. Close the "Settings" menu by tapping on the ⬆ button.
4. Start the recording of the video files [210].

6.7 Recording video files via SmartRecording

If SmartRecording is activated, video recordings of the live image are taken in in HD quality and stored in a buffer for five minutes. After five minutes have expired, the oldest video data in the buffer are cleared in such a way that the last five minutes always remain stored in the buffer. You can go back to the video files stored in the buffer and subsequently record them during an operation.

You can combine SmartRecording with video files in LowRes quality: Leave video recording in LowRes quality activated throughout the entire duration of use of the device. This recording requires much less storage space than a recording in HD quality. To record important surgical situations in HD quality, activate SmartRecording so that you can start the recording extending five minutes into the past for an important surgical situation at any time.

The buffer storage is automatically cleared and restarted if:

- The current user is changed
- The video recording is changed from 2D to 3D or vice versa
- The parallel recording is changed from "USB" to "Network" or vice versa

Activating SmartRecording

Action

1. Tap on  Settings  Video Recording.
2. Activate the "SmartRecording" slide switch.
 - ⇒ From this point on, the video files will be written to a buffer even if no video recording has been started.
 - ⇒ After five minutes have expired, the oldest video data in the buffer are cleared in such a way that the last five minutes always remain stored in the buffer.
3. Close the "Settings" menu by tapping on the  button.

Prerequisite

Starting SmartRecording

- ☑ No video streaming is running on the device
- ☑ No HD video recording is running on the device
- ☑ The internal hard drive has enough free space; the hard drive symbol in the status bar must not be colored orange

4. Tap on the [Video Recording] button  on the screen or actuate the preconfigured [Video Recording] button on the handgrip or foot control panel.

⇒ The runtime display in the "Video" field does not start at 00:00 as usual, but rather, at the current runtime of the buffer.

⇒ On the screen, the [Video Recording] button changes from  to .

⇒ A blue bar with information on the video recording appears at the bottom of the screen:

On the left side: HD recording; on the right side: storage location, e.g. internal hard drive with specification of storage capacity in percentage.

⇒ The blue bar with information on the video recording turns orange as soon as the free storage capacity of the internal hard drive is only 20%.

5. As soon as the blue bar with information on the video recording turns orange, stop the video recording and make sure that there is sufficient storage space on the internal hard drive.

6. To stop the video recording, tap on the [Video Recording] button  on the screen or actuate the preconfigured [Video Recording] button on the handgrip or foot control panel.

⇒ The video recording ends.

⇒ Recording of the video files to the buffer starts again.

Ending SmartRecording

7. Tap on  Settings  → Video Recording.

8. Deactivate the "SmartRecording" slide switch.

⇒ From this point on, video files will no longer be written to the buffer.

⇒ If you now start recording HD video files, the runtime begins at the moment you start the video recording (at runtime 00:00).

6.8 Streaming video files

The network can live stream video files via the following three network connections:

- LAN
- WLAN
- WLAN hotspot

You can use only one of the three options for streaming at any given time. As soon as the respective network connection has been stably established, the ZEISS KINEVO 900 sends the streaming signals to a streaming client via one of these three options. The streaming client then forwards the signal to the output device.

The streaming client must support the following protocols and video codecs:

- Protocol: RTSP (RTP via UDP)
- 2D video codec:
 - Encryption: H.264 / MPEG-4 AVC / High Profile / Level 4.2
 - HD resolution: 1920x1080p50/60
 - LowRes resolution: 960x540p25/30
- 3D video codec:
 - Encryption: H.264 / MPEG-4 AVC / High Profile / Level 4.2
 - Resolution: 1920x1080p50/60 TopBottom

If mobile systems are to be used for the display of the video stream, ZEISS recommends using the VLC media player (VideoLAN) for Apple iOS and Windows.

Prerequisite

- ☒ The desired network connection (LAN, WLAN, WLAN hotspot) has been set up by the IT admin and is stable
- ☒ The streaming client supports the above protocols and video codecs
- ☒ Data transmission rate > 40 Mbit/s

Action

1. Tap on  Settings  Video Recording.
2. Select the resolution for the video image: "HD Video" or "LowRes Quality".
3. Activate the slide switch in the "Video Streaming" field.
 - ⇒ As soon as video streaming has been activated, the address of the ZEISS KINEVO 900 appears in the "Web Address for Video Streaming" field.
4. Enter this address in the streaming client so that your output device can receive the streaming signal from the ZEISS KINEVO 900.
 - ⇒ As soon as the connection is stable, video streaming begins.

6.9 Taking photos with the integrated video camera

Photos are snapshots and can be taken in various situations:

- During an operation when no video recording is running
- During an operation when a video recording is running These photos are created as markers in the video files so that they can be located quickly during playback of the recorded video files.
- If, while playing back a recorded video file in the "Patient" menu, you notice an interesting location you would like to look at individually later on

The recorded photos are assigned to the current patient in each case.

Prerequisite

- ☒ The "Photo" function has been assigned to a button on the handgrip or the foot control panel

Action

1. Tap on  Settings →  Photo.
2. Select the "Internal Camera" function in the "Recording" field.
3. Select the desired format in the "Storage format (internal camera)" field.
4. Close the "Settings" menu by tapping on the  button.
5. Take a photo by tapping on the [Photo]  button on the screen or actuating the preconfigured [Photo] button on the handgrip or foot control panel.
 - ⇒ A thumbnail of the photo briefly appears at the bottom of the screen.
 - ⇒ An acoustic signal sounds.
 - ⇒ The photo is then saved in the "Patient" menu under the current patient.
6. Tap on the [Photos] button at the bottom of the screen or open the "Patient" menu to view the photos.

6.10 Taking photos with the DSLR camera (option)

You can optionally connect a Canon DSLR camera to the device via a photo adapter and take photos with it. Photos are snapshots and can be taken in various situations:

- During an operation when no video recording is running
- During an operation when a video recording is running

Note:

The photos taken are saved exclusively to the storage medium of the DSLR camera.

Prerequisite

- ☑ The microscope is configured as a system for hybrid (optical-digital) visualization.
- ☑ The photo adapter, T2 adapter and DSLR camera have been installed at the lateral co-observation port of the surgical microscope
- ☑ The DSLR camera and the connection socket for photo release are interconnected via the Canon EOS digital control cable (order no.: 000000-1229-877)
- ☑ The "Photo" function has been assigned to a button on the handgrip or the foot control panel

Action

1. Tap on ⚙ Settings → 📷 Photo.
2. Select the "External Camera" function in the "Recording" field.
3. Activate the maximum light output of 100% in the "Light intensity for external camera" field.
 - ⇒ The switch illuminates blue .
 - ⇒ The light intensity of the xenon illumination is briefly set to 100% when the external DSLR camera is released.
4. Close the "Settings" menu by tapping on the ⏮ button.
5. Take a photo by tapping on the [Photo] 📷 button on the screen or actuating the preconfigured [Photo] button on the handgrip or foot control panel.
 - ⇒ The photo is saved to the external DSLR camera.

6.11 Patient data, photographs and videos

CAUTION!

Do not use images and videos for diagnostic purposes!

The monitors are neither calibrated nor designed for diagnostic purposes.

- The video contents and images displayed on the monitors (cut sequences) must not be used for diagnostic purposes. They may contain deviations with respect to scale, shape, contrast and color.

6.11.1 Opening the patient directory

Action

1. Tap on the [Patient] button in the status bar.
 - ⇒ The patient directory opens. The data of the current patients are displayed.

6.11.2 Display of media objects (files)

Three views are available for displaying the files in the patient directory. Tapping on one of the symbols changes the number of files displayed.

6.11.3 Viewing patient images

Action

1. Tap on the [Photos] button in the "Patient" menu.
 - ⇒ All photos of the patient are displayed.
If not all of the photos are displayed, scroll down the photo list.
2. Tap on the photo you want to view in full size.
 - ⇒ The photo is displayed in full size.
The following information/buttons are displayed in the footer:
 - ⇒ Opened photo/number of photos
 - ⇒ Date and time of recording
 - ⇒ Located on the left and right edge of the photo are arrows [<], [>] for moving in reverse or forward.
3. Tap on the appropriate arrow [<], [>] to view the previous photo or the next photo or "wipe" horizontally across the touchscreen.
4. Close the menu by tapping on the [Patient] button or on  in the status bar.

Photos with position memory markings  are used to move to saved position memory positions.

6.11.4 Viewing patient videos

Action

1. Tap on the [Videos] button in the "Patient" menu.
 - ⇒ All recorded video files of the patient are displayed.
If not all of the videos are displayed, scroll up the video list.
2. Tap on the [Start] button  to start the playback. Tap on the [Pause] button  to stop the playback.
3. Tap on the [Fast Forward]  or [Fast Rewind]  button to run the video forward or rewind it at fast speed.
4. Tap on the [Photo]  button to capture and save a single image from the sequence. A marker is also set simultaneously to provide you with orientation in the video.
 - ⇒ An acoustic signal sounds.
 - ⇒ The photo is then saved in the "Patient" menu under the current patient.
 - ⇒ This location is furnished with a marker in the video file and marked with a small gray dot located above the playback bar in the video file.
5. To go to the locations where you created photos either while recording or while subsequently viewing the video file, move the "Video playback" slide control to any location where a small gray dot is displayed.


6. Tap on the [Edit video] button  to edit a recorded video [221].
7. Tap on the [Delete] button  to remove the current video.
8. Tap on the [Export] button  to export the current video.
9. Located on the left and right edge of the video are arrows [, ,] for moving in reverse or forward. Tap on the appropriate arrow [, ,] to open the previous video or the next video.

6.11.5 Editing patient videos

Videos must not be edited while a recording is in progress. Stop the recording first.

Action

1. Tap on the [Videos] button in the "Patient" menu.
2. Tap on the video to be edited.
⇒ The video view opens.
3. To edit the video, tap on the [Edit video] button .
⇒ The "Video Editing" menu opens.

4. Tap on the [Start]  button to play back the video recording.
5. Tap on the [Fast Forward]  button to run the video recording forward at fast speed. You also can move the slide control underneath the video recording to jump to a specific location.
6. When you have found the part of the video recording you would like to exclude from the summary, tap on the [Pause]  button.
⇒ The current playback position in the video recording is displayed in the timeline selection window as a time marker (small white bar).
7. Drag the blue start bar of the timeline selection window onto the time marker.
8. Tap on the [Cut in] button.
⇒ The video sequence starts at this position.
9. Tap on the [Start]  button to continue the playback of the video recording.
10. Tap on the [Pause] button  if you do not want to play back the entire video recording, but would like to remove the rear portion of the video from the display.
⇒ The current playback position in the video recording is displayed as a time marker (small white bar) in the timeline selection window.
11. Drag the blue end bar of the timeline selection window onto the time marker.

12. Tap on the [Cut out] button.
⇒ The video sequence now ends earlier.
13. You can relocate the start and end points of the video sequence at any time.
14. Tap on the [Temporarily save video sequence] **[+]** button to temporarily save the marked video sequence.
⇒ The cut-out video footage is displayed as a preview at the bottom of the screen.
15. Tap on the [Create New Video] button to save it in the "Patient" menu.
16. To load and edit another video of the patient, tap on the [Add original recording] + button.
⇒ A preview window with all original video recordings of the patient opens.
17. Tap on the video recording you want to edit.
18. Edit the video recording as described.
19. Close the menu by tapping on the [Patient] button or on  in the status bar.

6.12 DICOM (option)

Transfer of patient data

The DICOM option enables a standardized exchange of patient data with the PACS system of the clinic.

Please observe the DICOM Conformance Statement applicable to the KINEVO 900: G-30-1952.

In order to utilize the DICOM functionality, the KINEVO 900 must be registered and configured in the hospital network.

NOTE

Do not use imported images and videos for diagnostic purposes!

Perform diagnosis only at approved diagnostic stations.

Playback images may contain deviations with respect to scale, shape, contrast and color.

- ▶ It is best not to enter patient master data locally in order to avoid mix-ups.
- ▶ Preferably load patient master data from a DICOM Worklist.
- ▶ If data are transmitted by the device into the data network or are kept available in the data network, there is a risk that they will be corrupted or transmitted incompletely. Therefore, no liability can be accepted for the correctness of the data.

6.12.1 Importing patient data from the DICOM Worklist

Action

1. Tap on the [Patient] button in the status bar.
⇒ The "Patient" menu is displayed.
2. Tap on the [Import]  button.
3. Tap on the [DICOM] button.
4. Fill out the offered search criteria.
The search criteria are interpreted as "AND conjunction". Empty search windows are interpreted as "wild cards" (placeholders for other characters).
To limit the search, you can e.g. enter the first letters of the patient's name.
5. Select the scheduled date.
6. Select the [Modality].
7. Tap on the [Display worklist] button.
⇒ The DICOM worklist with the planned patient data records is shown.
8. Accept the patient record by marking the box in front of the patient name  or all records by marking the [planned for]  box.
9. Tap on the [Import Selected] button in order to import the data records.

6.12.2 Exporting patient data to DICOM

Action

1. Tap on the [Patient] button in the status bar.
⇒ The "Patient" menu is displayed.
2. Mark all data to be imported .
3. Tap on the [Export Selected]  button at the bottom right on the touchscreen.
⇒ The "Export Patient Data" menu is displayed.
4. Tap on the [DICOM] selection field in the "Export Path" field.
5. Tap on the [Start] button to export the patient data.
⇒ The message "Export in progress" is displayed.

6.13 Use QEVO Micro-Inspection Tool (option)

The Micro-Inspection Tool QEVO is intended for the inspection of open surgical fields in general surgical interventions and for the visual imaging of ventricles and structures inside the brain during neurological surgical interventions.

When connecting and using the Micro-Inspection Tool, is vital that you observe the relevant warning notices and safety information in the QEVO and QEVO ECU Instructions for Use! [13]

6.13.1 Switch QEVO on / off

Prerequisite

- ☒ ZEISS KINEVO 900 is switched on and ready for use.
- ☒ QEVO is connected to ZEISS KINEVO 900 . [167]
- ☒ QEVO is activated on ZEISS KINEVO 900 . [168]
- ☒ A button on the hand grip and/or on the foot control panel of the ZEISS KINEVO 900 is configured for switching on QEVO . [180]

Action

- ▶ To switch on QEVO press the [QEVO] button on the hand grip or on the foot control panel.
 - ⇒ Light und camera of the QEVO switch themselves on.
 - ⇒ The camera image of QEVO is shown as a full image. The camera image of ZEISS KINEVO 900 is displayed on the touchscreen as Picture in Picture.
- ▶ Press the [QEVO] button on the hand grip or on the foot control panel to switch off QEVO.
 - ⇒ Light and camera of the QEVO switch off.
 - ⇒ The camera image of ZEISS KINEVO 900 is shown.


6.13.2 Display QEVO user interface

Prerequisite

☒ QEVO is switched on. [224]


Action

1. Tap on the touchscreen.
⇒ The control elements of the ZEISS KINEVO 900 software are shown.
2. Tap on the bar on the right-hand edge of the screen.

⇒ Additional menu "Extended operation" opens.

3. Tap on the [QEVO] button in the upper area of the additional menu.

⇒ The "QEVO" menu opens.

6.13.3 Switching the automatic switch-off function on/off

The switch-off function switches QEVO off if the tool is set aside and not moved for approx. 40 seconds. If the automatic switch-off function is not switched on, the light source remains switched on continuously. QEVO must then be switched off by pressing the [QEVO] button on the hand grip or the foot control panel.

Prerequisite

- ☑ The "QEVO" menu [225] or the "MultiVision" menu [106] is opened.

Action

- ▶ To switch on the automatic switch-off function, activate the [Automatic Switch-off] button.
 - ⇒ The button will light up in blue .
- ▶ To switch off the automatic switch-off function, deactivate the [Automatic Switch-off] button.
 - ⇒ The button will turn gray .


6.13.4 Perform white balance

Prerequisite

- ☑ QEVO is switched on.
- ☑ The "QEVO" menu is opened.
- ☑ A white surface is shown.

Action

1. Hold the shaft tip at a distance of 5 to 30 mm on a white surface.

2. To start the white balance, tap on the [White balance] button .

6.13.5 Setting the image brightness

QEVO regulates the image brightness automatically. If the shaft tip is directed on a light surface, the image brightness is reduced. If the shaft tip is then introduced into the surgical field quickly, it may take 1 or 2 seconds until you see a fully illuminated image.


The image brightness can be set in two different ways.

Setting the image brightness on the hand grip

In the QEVOmode, both yellow buttons of the right hand grip are automatically reserved for controlling the image brightness of the Micro-Inspection Tool.

Prerequisite

- ☒ QEVO is switched on.

1	[Increase brightness] button	2	[Reduce brightness] button
---	------------------------------	---	----------------------------


Action

1. To increase the image brightness, press the yellow button on the right-hand side of the right hand grip often enough until the desired image brightness is set or the maximum image brightness is reached.
2. To reduce the image brightness, press the yellow button on the left-hand side of the right hand grip often enough until the desired image brightness is set or the minimum image brightness is reached.

Prerequisite

Setting the image brightness in the "QEVO" menu

- ☒ QEVO is switched on.
- ☒ The "QEVO" menu is opened. [225]
- 3. In the "QEVO" menu, tap on [Brightness].
⇒ A slide control opens.

- 4. To increase the image brightness, push the slide control to the right or tap $\oplus$ often enough that the desired image brightness is set or the maximum image brightness is reached.
- 5. To decrease the image brightness, push the slide control to the left or tap $\ominus$ often enough that the desired image brightness is set or the minimum image brightness is reached.

6.13.6 Change QEVO alignment

Using the "QEVO alignment" function, you can rotate the QEVOcamera image on the monitor in order to align the viewing direction of the Micro-Inspection Tool relative to the viewing direction of the surgical microscope.

The alignment can be set in increments of 45°. The blue segment in the "QEVO" menu symbolizes the current alignment of the instrument relative to the viewing direction of the surgical microscope.

The QEVO alignment can be set in two different ways.

Prerequisite

Set QEVO alignment on the handgrip

- ☒ QEVO is switched on.


1	XY Joystick
---	-------------

Action

- In order to align the image in a clockwise direction, press to the right on the XY joystick of the right hand grip often enough until the image corresponds to the alignment of the QEVO.
- In order to align the image in a counterclockwise direction, press to the left on the XY joystick of the right hand grip often enough until the image corresponds to the alignment of the QEVO.

*Prerequisite***Set QEVO alignment in the "QEVO" menu**


- ☒ QEVO is switched on.
- ☒ The "QEVO" menu is opened. [225]

- In order to align the image in a clockwise direction, tap  often enough until the image corresponds to the QEVO alignment. Alternatively, tap on the corresponding circular area of the section.
- In order to align the image in a counterclockwise direction, tap  often enough until the image corresponds to the QEVO alignment. Alternatively, tap on the corresponding circular area of the section.
- To reset the image rotation, tap on . Alternatively, tap on the upper section of the circular area.

6.13.7 Show / hide marker

The marker shows the upper edge of the recorded camera image. The position of the marker changes accordingly with the rotation of the Micro-Inspection Tool around the optical axis.

1	shown marker
---	--------------

Prerequisite

- ☒ QEVO is switched on.
- ☒ The "QEVO" menu [225] or the "MultiVision" menu [106] is opened.

Action

- In order to show the marker, activate the [Marker] button.
⇒ The button will light up in blue .
- In order to hide the marker, deactivate the [Marker] button.
⇒ The button will turn gray .

6.13.8 Switching the light on/off

If needed, you can use QEVO as an additional light source for the surgical field, by switching on the light of the Micro-Inspection Tool without camera. To do this, a button on the hand grips or on the foot control panel must be configured for switching the QEVO light on and off.

Prerequisite

- ☒ A button on the hand grip and / or on the foot control panel of the ZEISS KINEVO 900 is configured for switching on the QEVO light . [180]

Action

- To switch QEVO light source on, press the [QEVO light ON/OFF] button on the hand grip or on the foot control panel.
⇒ The QEVO light switches on.
- To switch off the QEVO light, press the [QEVO light ON/OFF] button on the hand grip or on the foot control panel.
⇒ The light of the QEVO switches off.

6.14 BLUE 400 (option)

NOTE

Application of BLUE 400 only with eyepiece-based optical observation!

Tubes and eyepieces must be used when the BLUE 400 option is applied.

Note:

Note that if the user is color blind, this may impact the application of BLUE 400 (recognizing and distinguishing red tones in fluorescent tissue).

Because of regulatory restrictions regarding the contrast medium used with BLUE 400, the BLUE 400 option is not registered in all countries, and is therefore not available in every country. Please contact your local ZEISS representative to obtain information regarding the registration status of BLUE 400 in your country. In the USA, the BLUE 400 option may not be used for clinical applications, but only for research purposes!

Intended use / Indication for use (US only)

The BLUE 400 option is designed for a surgical microscope and can be used in tumor surgery to visualize marked tumor tissue with fluorescence dyes.

Normal use

The integrated fluorescence module BLUE 400 serves to render fluorescent areas visible. It is designed for excitation in the wavelength range of 400 to 410 nm and observation in the wavelength range of 620 to 710 nm. The surgeon can switch back and forth between white light and blue excitation light for fluorescence at the touch of a button. For optimal viewing of the fluorescent light, the room or the ambient lighting should be dimmed.

Normal use under standard conditions

The BLUE 400 fluorescence option enables the user to visualize fluorescence light and auto fluorescence light from tissue and to display it digitally.

The visualization of the fluorescence signal is influenced by several different factors:

- the fluorescence medium and its concentration in the tissue
- the illumination intensity of the light source in the defined wavelength range
- the transmission of the optical system
- the total magnification and the aperture setting
- the working distance and the illuminated field diameter

When using the KINEVO 900 in conjunction with a fluorescence substance, the medical contraindications apply as in the use of suitable marker substances.

- Only use fluorescent agents approved for the planned application.
- Check whether the fluorescent agent can be excited in the wavelength range from 400 to 410 nm and whether it emits fluorescent light of sufficient intensity in the wavelength range from 620 to 710 nm.
- As in almost all procedures, false-positive and false-negative results can also occur in the fluorescence-based method. Evaluation by the user based on other methods is necessary.

For optimal visualization of the fluorescence signal, the room and surrounding lights should be dimmed.

6.14.1 Configuring BLUE 400

6.14.1.1 Configure fluorescence button on the hand grip or FCP

Action

1. Tap on  Settings →  Operation → Hand grips [180] or FCP → Configure fluorescence button.
2. Tap on a configurable button.
 - ⇒ The available functions appear in the selection field on the right side.
3. Tap on the "Fluorescence" function in the selection field.
 - ⇒ The function is accepted and displayed next to the button symbol.

6.14.1.2 Configuring the light intensity for BLUE 400

Action

1. Tap on  Settings →  Fluorescence → BLUE 400.
2. Set the desired start value with the [Light Intensity (start value)] slide control in the "BLUE 400" field.
 - ⇒ The illumination lights with this brightness when BLUE 400 is switched on.

In order to achieve a good fluorescence result, ZEISS recommends using the factory setting of 100%.

6.14.1.3 Configuring the video camera for BLUE 400

Action

1. Tap on  Settings →  Video Camera.
2. Scroll down in the Video Camera menu to → BLUE 400.
3. Set the shutter control to "Auto" or "Manual".
4. If you have set the shutter control to "Auto", set the desired brightness level on the [Auto Brightness] controller.
5. If you have set the shutter control to "Manual", set the exposure time in the [Manual Shutter] selection field.
6. Close the "Settings" menu by tapping on the .

6.14.2 Carry out BLUE 400 functional inspection

CAUTION!

Risk of infection due to unsterile fluorescence target.

The patient may become infected if an unsterile fluorescence target is used.

- ▶ Therefore, perform the fluorescence function test prior to the operation and without a patient.

CAUTION!

Do not use images and videos for diagnostic purposes!

The monitors are neither calibrated nor designed for diagnostic purposes.

- ▶ The video contents and images displayed on the monitors (cut sequences) must not be used for diagnostic purposes. They may contain deviations with respect to scale, shape, contrast and color.

Prerequisite

- ☒ BLUE 400 fluorescence target is available

Action

1. Place the BLUE 400 fluorescence target in the center of the illuminated field.
2. Set the focus to 250 mm.
3. Set the zoom to 3 x.
4. Focus the fluorescence target by moving the surgical microscope up and down.
5. Dim the room lighting.
6. Start the BLUE 400 fluorescence application [235].
7. Check whether the image on the touchscreen corresponds to the supplied sample image in the target.
NOTE! There may be some slight differences in color and brightness.
8. Stop the BLUE 400 fluorescence application.

6.14.3 Starting the BLUE 400 fluorescence application

As soon as the BLUE 400 fluorescence application is started, the "MultiVision" function is automatically deactivated in order to prevent any impairment of the image quality in the eyepiece and on the video image resulting from the additional data injection.

Prerequisite

- ☑ In the "Fluorescence" menu, BLUE 400 is selected in the "Fluorescence mode" field.

Action

1. To start the BLUE 400 fluorescence application, tap on the handgrip or the FCP button to which you have assigned the "Fluorescence" function.
2. Or, in the main menu bar on the touchscreen, tap on  Fluorescence → Select BLUE 400 → Start the fluorescence application → .

Result

- ✓ The BLUE 400 fluorescence application begins.

6.15 YELLOW 560 (option)

NOTE

Application of YELLOW 560 only with eyepiece-based optical observation!

Tubes and eyepieces must be used when the YELLOW 560 option is applied.

Because of regulatory restrictions regarding the contrast medium used with YELLOW 560, the YELLOW 560 option is not registered in all countries, and is therefore not available in every country. Please contact your local ZEISS representative to obtain information regarding the registration status of YELLOW 560 in your country.

Intended use / Indication for use (US only)

The YELLOW 560 option is an accessory for the intraoperative illumination of the surgical area with light in the wavelength range from 460 nm to 500 nm. Light in the wavelength range between 540 and 690 nm is particularly highlighted.

Normal use

The integrated YELLOW 560 fluorescence module is used to make fluorescing areas visible. It is designed for excitation in the wavelength range of 460 to 500 nm and observation in the wavelength range of 540 to 690 nm. The surgeon can switch back and forth between white light and blue excitation light for fluores-

cence at the touch of a button. For optimal viewing of the fluorescent light, the room or the ambient lighting should be dimmed.

Normal use under standard conditions

The YELLOW 560 fluorescence option permits the user to visualize and digitally record fluorescent light and autofluorescent light emitted by tissue.

The visualization of the fluorescence signal is influenced by several different factors:

- the fluorescence medium and its concentration in the tissue
- the illumination intensity of the light source in the defined wavelength range
- the transmission of the optical system
- the total magnification and the aperture setting
- the working distance and the illuminated field diameter

When using the KINEVO 900 in conjunction with a fluorescence substance, the medical contraindications apply as in the use of suitable marker substances.

In almost all procedures, even the fluorescence-supported procedures, false-positive and false-negative results can occur.

For optimal visualization of the fluorescence signal, the room and surrounding lights should be dimmed.

6.15.1 Configuring YELLOW 560

6.15.1.1 Configure fluorescence button on the hand grip or FCP

Action

1. Tap on  Settings →  Operation → Hand grips or FCP → Configure fluorescence button.
2. Tap on a configurable button [180].
 - ⇒ The available functions appear in the selection field on the right side.
3. Tap on the "Fluorescence" function in the selection field.
 - ⇒ The function is accepted and displayed next to the button symbol.

6.15.1.2 Configuring the light intensity for YELLOW 560

In order to achieve a good fluorescence result, ZEISS recommends using the factory setting of 100%.

Action

1. Tap on  Settings →  Fluorescence → YELLOW 560.
2. Set the desired start value with the YELLOW 560 [Light Intensity] controller.
 - ⇒ The surgical field will be illuminated with this light intensity as soon as YELLOW 560 is activated.

* In China the light intensity cannot be set to 100% in the YELLOW 560 fluorescence mode due to national regulations.

6.15.1.3 Configuring Video Recording for YELLOW 560

Action

1. Tap on  Settings →  Video Camera.
2. Scroll down in the "Video Camera" menu to → YELLOW 560.
3. Set the shutter control to "Auto" or "Manual".
4. If you have set the shutter control to "Auto", set the desired brightness level on the [Auto Brightness] slide control.
5. If you have set the shutter control to "Manual", set the desired exposure time in the [Manual Shutter] selection field.
6. Close the "Settings" menu by tapping on the  button.

6.15.2 Procedure for YELLOW 560 function test

CAUTION!

Risk of infection due to unsterile fluorescence target.

The patient may become infected if an unsterile fluorescence target is used.

- Therefore, perform the fluorescence function test prior to the operation and without a patient.

CAUTION!

Do not use images and videos for diagnostic purposes!

The monitors are neither calibrated nor designed for diagnostic purposes.

- The video contents and images displayed on the monitors (cut sequences) must not be used for diagnostic purposes. They may contain deviations with respect to scale, shape, contrast and color.

Prerequisite

- ☒ YELLOW 560 fluorescence target is available

Action

1. Place the YELLOW 560 fluorescence target in the center of the illuminated field.
2. Set the focus to 300 mm.
3. Set the zoom to 4.5x.

4. Focus the fluorescence target by moving the surgical microscope up and down.
5. Dim the room lighting.
6. Start the YELLOW 560 fluorescence application [▶ 238].
7. Check whether the image on the touchscreen corresponds to the supplied sample image in the target.
NOTE! There may be some slight differences in color and brightness.
8. Stop the YELLOW 560 fluorescence application.

6.15.3 Starting the YELLOW 560 fluorescence application

As soon as the YELLOW 560 fluorescence application is started, the "MultiVision" function is automatically deactivated in order to prevent any impairment of the image quality in the eyepiece and on the video image resulting from the additional data injection.

Prerequisite

- ☑ In the "Fluorescence" menu, YELLOW 560 is selected in the "Fluorescence mode" field.

Action

1. To start the YELLOW 560 fluorescence application, tap on the handgrip or the FCP button to which you have assigned the "Fluorescence" function.
2. Or tap on fluorescence → YELLOW 560 → Starting the fluorescence application →  in the main menu list on the touchscreen.

Result

- ✓ The YELLOW 560 fluorescence application starts.

6.16 Using the device as a hotspot

As soon as the WLAN adapter has been activated, you can use the device as a hotspot.

After each reboot of the device, the hotspot must be switched on again by the "Default Users" or "Users" user group. When the device is shut down, the hotspot function also is switched off.

The WLAN option required to use the device as a hotspot is not available in all countries.

Prerequisite

- ☑ The device is connected to the network via WLAN [151].
- ☑ Device is created as a hotspot and activated [151].


Action

1. Tap on the [Status Information] button in the status bar.

⇒ The "Status Information" menu is displayed.

2. Activate the [Hotspot] slide switch.

Result

- ✓ The device is available as a hotspot.

6.17 Connecting an external device to the device via web interface (option)

If the web interface is activated in the device, you can access the patient data from an external device.


The WLAN option required to use the device as a hotspot is not available in all countries.

Prerequisite

- ☒ The device is created as a hotspot and activated [151].
- ☒ Device is being used as a hotspot [239].
- ☒ Web interface of the the device is activated [152].

Action

1. Tap on the [Status Information] button in the status bar.
 - ⇒ The "Status Information" menu is displayed.
 - ⇒ The address of the web interface is visible.

2. Enter the displayed address of the web interface on the external device.
 - ⇒ An input mask for the web interface password opens on the external device.
3. Enter the web interface password (assigned by the "IT Admin" user group) on the external device.
 - ✓ The external device is connected to the device via the web interface; patient data can be accessed.

Result

6.18 Using ZEISS Smart Services

The following settings and configurations may only be carried out by authorized users with IT admin system rights [149].

6.18.1 Activating ZEISS Smart Services

Prerequisite

- ☑ You are logged on as Administrator.

Action

1. Tap on ⚙ Settings → Connections.
2. Scroll to the paragraph titled "ZEISS Smart Services".
3. Activate the slide switch in the "Activate" field.
⇒ In the "Extras" [109] menu, the "Diagnostic Package" and "Remote Service Access" functions become available.

6.18.2 Send Diagnostic Package

Prerequisite

- ☑ The device is connected to the network via LAN or WLAN. [149]

Action

- ☑ The ZEISS Smart Services have been activated. [241]
1. Tap on ⚙ Settings → Extras.
 2. Tap on the [Send Diagnostic Package] button.
⇒ A data transfer progress bar is displayed on the monitor.

6.18.3 Allow Remote Service Access

CAUTION!

Risk of injury due to accidental device travel movements during remote maintenance!

If ZEISS Smart Services is activated, the device can be controlled and possibly moved via remote maintenance.

- ▶ Make sure that the network connection is stable.
- ▶ Make sure that the device is not standing within reach of a patient.
- ▶ Make sure that the device is not standing inside a sterile surgical area.
- ▶ Do not leave the device unattended when ZEISS Smart Services is activated.
- ▶ Reboot the device before using it for clinical applications.

Prerequisite

- ☑ The device is connected to the network via LAN or WLAN. [149]

Action

- ☑ The ZEISS Smart Services have been activated. [241]
1. Tap on ⚙ Settings → Extras.
 2. Activate the [Remote Service Access] slide switch.
 3. Deactivate the [Remote Service Access] slide switch as soon as the remote service has been completed.

6.19 Powering the device down

Action

1. Press the Standby/ON-OFF switch on the device.
 - ⇒ The device shuts down
 - ⇒ The device is in Standby mode if the Standby/ON-OFF switch illuminates white.
2. Remove the power plug from the power outlet to safely disconnect the device from the power supply.

7 Cleaning and disinfection

7.1 Safety during cleaning and disinfection

WARNING!

Contamination with pathogenic germs!

If the device is cleaned or disinfected using unsuitable methods or cleaning agents, pathogenic germs can accumulate on the device and the device may also be damaged.

- ▶ Make sure that the device is operated only in a sterile manner and with sterile accessories, Attaching a ZEISS SMARTDRAPE.
- ▶ Dispose of sterile accessories in line with the applicable regulations for hospital wastes and/or materials contaminated by bacterial and/or viral infections.
- ▶ Read and follow the Instructions for Use for the accessories used.
- ▶ Make sure that the device is operated, cleaned and disinfected only by personnel who have been instructed in these procedures.
- ▶ When cleaning and disinfecting the device, only use cleaning/disinfection agents that are suitable and approved for the surface in question.

7.2 Contamination of the system

Dust can penetrate into the internal optics of the device or its individual components.

Action

- ▶ Always close all openings that are not in use with the covers provided (e.g. openings for eyepieces, tube outlets or the lateral co-observation outlets).
- ▶ Always store tubes, eyepieces and accessories in dust-free cases when they are not being used.
- ▶ After use, cover the system with the supplied dust cover to protect it from dust.
- ▶ Clean used accessories directly after use.

7.3 Cleaning and disinfecting agent

7.3.1 Cleaning agents

	Unit	Order number
Optics cleaning kit	1 piece	000000-2096-685
Microfiber cleaning cloth	1 piece	000000-1439-266
Optical lens wipes	Available from specialized dealers	

7.4 Cleaning

Clean the device before the first use and after every use.

7.4.1 Cleaning optical surfaces

The multi-layer T* coating of the optical components (e.g. eyepieces, objective lenses) ensures optimal image quality. Image quality is impaired by even slight contamination of the optics or by a fingerprint. Clean the exterior surfaces of the optical components (eyepieces, objective lenses) only when necessary:

- ▶ Do not use any chemical cleaning agents.
- ▶ Use a clean and grease-free brush to remove dust.

TIP: For the regular cleaning of surgical microscope's objective lenses and eyepieces, we recommend the optics cleaning set and Cleaning agents [244] available from ZEISS.

7.4.2 Cleaning the monitor (touchscreen)

To prevent damage, ensure that no moisture or cleaning agents get inside the monitor.

NOTE

Incorrect use of glass cleaner!

Damage to the monitor.

- ▶ Never spray glass cleaner directly onto the monitor.

☒ Switch off the device.

- ▶ Clean the monitor with a soft, clean microfiber cloth or with moist optical lens wipes.
- ▶ If necessary, moisten the cloth slightly with water or pure glass cleaner.
- ▶ Wipe the monitor with a moist cloth.

Action

Prerequisite

Action

7.4.3 Cleaning mechanical surfaces

All mechanical surfaces of the system can be cleaned by wiping them with a damp cloth.

Action

- ▶ Do not use any aggressive or abrasive cleaning agents.
- ▶ Remove any possible residue using a mixture of 50% ethyl alcohol and 50% distilled water plus a dash of household dishwashing liquid.

7.4.4 Fogging of optical surfaces

We recommend using an anti-fogging agent to prevent fogging of optical surfaces. Anti-fogging agents like the ones offered by opticians for applications with eyeglasses are also suitable for optical surfaces from ZEISS.

Action

- ▶ Observe the Instructions for Use pertaining to the anti-fogging agent concerned.

An anti-fogging agent does not ensure fog-free eyepiece optics. It cleans eyepiece optics and protects them against dirt, grease, dust, lint and fingerprints.

7.5 Disinfection

7.5.1 Disinfecting the mechanical surfaces

The maximum concentrations are:

- For alcohol (tested with isopropyl alcohol): 60%
- For aldehyde (tested with glutaraldehyde): 2%
- For quaternary compounds (tested with DDAC): 0.2 %

To achieve maximum disinfection, disinfectants with an alcohol concentration higher than 70% can be used. There is a chance that long-term use of such disinfectants may cause the device's surfaces to become dull or matte, or may loosen the adhesive labels on the device without causing them to fall off. However, the use of such disinfectants will never impair device performance or endanger patients.

NOTE

Surface damage caused by wrong disinfectants!

Performing disinfection with the wrong disinfectants may result in damage to the surfaces of the device.

- ▶ Use an aldehyde and/or alcohol-based disinfectant. The addition of quaternary compounds is acceptable.
- ▶ In order to prevent surface tensions, you may use only the disinfecting components specified above.

Action

- ▶ Disinfect all of the required surfaces.

Empty page, for your notes

8 Maintenance

8.1 Safety during maintenance

Regularly checking and maintaining the device is necessary to ensure safe and correct operation of the system and a long service life. For example, within the framework of a service and maintenance contract.

Device-specific technical knowledge is required for all maintenance work. It is frequently necessary to open the system for this work. Please contact your local ZEISS Service organization in due time for the performance of these maintenance activities.

You can find the ZEISS contact partner for your country on the Internet at the following website: www.zeiss.com/med

8.2 Maintenance schedule for the operator

8.2.1 User maintenance intervals

Component	Activity or test criterion
Electromagnetic compatibility (EMC)	► No regular inspections or maintenance are required in order to maintain electromagnetic compatibility (EMC).
According to country-specific time specifications	
Safety inspection	► Checking the device and electrical safety
Prior to each use	
Brakes	► Electrical and mechanical function test: Release brakes completely; braking force holds microscope securely in balanced state without any change of position.
Lamps	► Check: ■ Operating hours within permissible range ■ Replacement lamp available and ready for use We recommend replacing the lamp after: The specified operating hours have been reached (500 hours)

8.3 Maintenance schedule for the authorized service

8.3.1 Service maintenance intervals

Time Interval	Activity or test criterion
Every 12 months:	Basic device functions Specific function checks: Microscope, illumination, stand system, visual inspection of device identification, cleaning and inspection of optics, video functionality
Every 4 years:	Extended device functions Extended function checks: Microscope, illumination, stand system, visual inspection of device identification and device safety, cleaning and inspection of optics, video functionality, equipment-specific functional tests and checks of various options

8.3.2 Performance of safety inspections

The safety inspection (Technical safety test) serves to determine and evaluate product safety. The user of this product is obligated to perform and document safety inspections in compliance with IEC 62353.

CAUTION!

Risk of injury without safety inspection (technical safety test)!

Hazards and product deficiencies will not be detected in time and can have a negative effect on patients, users, or others.

- ▶ Have a safety inspection in compliance with IEC 62353 performed on time and to the prescribed extent. At the same time, be sure to comply with the relevant national regulations.

Prerequisite

- ☒ Only the manufacturer or qualified persons may perform the safety inspection.

Action

- ▶ Check to make sure the Instructions for Use are available.
- ▶ Check to make sure all marks and inscriptions on the device are legible.
- ▶ Check the leakage current and the protective ground conductors.
- ▶ Check to make sure the casters and locking devices function properly and are not worn-out.
- ▶ Perform a function check of all switches, buttons, sockets and indicator lamps of the system.

8.4 Remote service

ZEISS Smart Services offers a direct data exchange between the device and ZEISS Service for remote maintenance or troubleshooting.

For informationen on the activation and use of ZEISS Smart Services, read Using ZEISS Smart Services [241].

Empty page, for your notes

9 Troubleshooting

9.1 Localize faults

If a particular event occurs during operation or if the system registers an irregularity, the system shows you this disruption via a notification on the monitor.

Faults are saved in a service file. You can export this log file and transmit it to your ZEISS Service organization.

9.1.1 React to disruptions with notifications

Error notifications provide information on:

- Operational steps that are currently being carried out, in which the displayed error is occurring
- Suggested solution to remedy the displayed error
- If appropriate, a hint that the log files must be exported and sent to ZEISS Service

Action

1. If an error notification appears on the monitor, read the notification text carefully.
2. Tap the "Close" button or press the center of the joystick button on the left hand grip to close the error message.
3. Remedy the displayed error.
4. If this error message remains displayed on the monitor, export the log file and send it to ZEISS Service [251].

9.1.2 Export Log Files

Disruptions with notifications are saved in log files on the system. You can export these log files to a USB storage device and send them to ZEISS Service for error analysis.

Action

1. Connect a USB storage device to the USB mouth switch socket.
2. Tap on  Settings →  Extras → Export log files → Start.
⇒ The log files are then exported to the connected USB storage device.
3. Tap on the "Status Information" field in the status bar.
4. Tap on the [Eject] button in the "USB" field.
5. Remove the USB storage device.
6. Send these log files to your ZEISS Service organization.

9.1.3 Service information

You can find the ZEISS contact partner for your country on the following website: www.zeiss.com/med

9.2 Remote service

ZEISS Smart Services offers a direct data exchange between the device and ZEISS Service for remote maintenance or troubleshooting.

For information on the activation and use of ZEISS Smart Services, read Using ZEISS Smart Services [241].

9.3 Malfunctions without notifications

Fault	Cause	Remedy
No function at all	The device is not connected to the power supply.	► Connect the device to the power supply.
	The device is not switched on.	► Press the "Device Power On/Off" operating button once. The operating button is illuminated white.
	Automatic circuit breaker in power switch of stand responds.	► Press the "Device Power On/Off" operating button again.
No surgical field illumination on the microscope	Lamps 1 and 2 have failed.	► Replace both lamp containers [257].
	Failure of device electronics.	► Illuminate the OR field with an additional OR lamp. ► Inform ZEISS Service.
	Light source has not been switched on.	► Activate the "Light" function in the main menu on the monitor.
The surgical field illumination is too dark.	The set brightness is too low.	► Increase the brightness using the preconfigured buttons on the hand grip / foot control panel / rocker foot switch.
	The lamp is too old / weak.	► Replace the lamp container [257].
The surgical field is not evenly illuminated.	The light guide is defective.	► Have the light guide replaced by ZEISS Service.
The surgical field illumination is too bright.	The set brightness is too high.	► Reduce the brightness using the preconfigured buttons on the hand grip / foot control panel / rocker foot switch.
The focus cannot be set via the hand grip / foot control panel / rocker foot switch.	The system control is defective.	► Adjust the focus manually using the corresponding rotary knob on the microscope body [59].
The zoom cannot be set via the hand grip / foot control panel / rocker foot switch.	The system control is defective.	► Adjust the zoom manually using the corresponding rotary knob on the microscope body [59].

Fault	Cause	Remedy
Poor image quality	The drape cover glass on the objective lens is dirty or of poor optical quality.	► Clean the cover glass. ► Use the ZEISS SMARTDRAPE.
The monitor is black.	The system control is defective.	► Do not touch the monitor! This is necessary to ensure that you do not execute any unwanted actions. ► Switch the device off, wait briefly (for approx. 2 min) and then switch it back on. ► If the monitor is still black, contact ZEISS Service.
No video image appears on the external monitor	The external monitor is not connected to the device.	► Connect the external monitor to the corresponding video output of the device [148].
External video source is not displayed	The external video source is not connected to the device.	► Connect the external video source to the corresponding video input of the device.
The co-observer tube cannot be repositioned	The system control is defective.	► Adjust the pivoting mirror for the co-observer tube manually using the corresponding rotary knob on the microscope body [59].
The brakes of the axes are closed.	The line voltage has failed, the device is no longer being supplied with power.	► Pull and push [255] the microscope body and stand arms into the desired position manually.
Malfunction / failure of wireless foot control panel	Batteries are depleted.	► Replace the batteries in the foot control panel; read and follow the Instructions for Use supplied for the foot control panel, G-30-1706.
	The foot control panel is not completely paired with the device.	► Pair the foot control panel with the device again [129].

9.4 Basic Function Mode following system control failure

If the system control fails and the line voltage is still present, all filters in the light source and in the microscope will be swiveled out. The Focus Light Link and Auto Brightness functions are not available.

The button assignments of the hand grip, foot control panel and rocker foot switch last set remain stored.

WARNING!

Failure of digital image

If the 3D video system fails and the image is no longer displayed on the 3D monitor, the 3D video system may no longer be used for the operation.

- ▶ Complete the operation with a backup device or with an alternative means of visualization (fully digital configuration: loupe; hybrid configuration: tube).
- ▶ Unplug the device from the electric power supply and contact ZEISS Service.

CAUTION!

Limited illumination functions!

Following failure of the system control, the light intensity is decreased to approx. 50%. The "Focus Light Link" and "Auto Brightness" functions are not available.

- ▶ You can once again increase or further decrease the light intensity at any time in order to complete your surgical procedure.

You can set the following functions in the Basic Function Mode:

Function	Adjustments following failure of the system control
Release all axes	▶ Press the bottom button [AB] on the back of the hand grips.
Release the stand axes	▶ Press the top button [SB] on the back of the hand grips.
Focus +/-	▶ Adjust the focus using the buttons configured for this purpose on the hand grip / foot control panel / rocker foot switch.
Zoom +/-	▶ Adjust the zoom using the buttons configured for this purpose on the hand grip / foot control panel / rocker foot switch.
Light +/-	▶ Adjust the brightness using the buttons configured for this purpose on the hand grip / foot control panel / rocker foot switch.

9.5 Device in the de-energized state

If the line voltage fails and the device no longer is being supplied with power, the brakes of all axes close. You can move the individual axes of the stand manually by overcoming the braking effect.

Action


- ▶ To move the microscope, grasp the microscope body (and not the handgrips!) with both hands and pull or push the microscope body into the desired position.
- ▶ To execute movement in individual stand axes, push or press the vertical / horizontal arm and/or the microscope suspension into the desired position.

9.6 Automatic lamp change

If the device detects a defective lamp or a lamp fails during operation, an automatic lamp change takes place. You also can change the lamp manually on the monitor.

Action

1. Tap on ⚙ Settings → Microscope → Light.
2. In the "Remaining lamp life" paragraph, tap on the option button for the lamp to be activated.

⇒ The lamp is activated.

3. Close the "Settings" menu by tapping on the ⌂ button.
⇒ The live image of the video camera appears on the monitor.
4. The light intensity of new lamps is higher; readjust the light intensity if necessary.

9.7 Manual lamp change

In case of motor failure


Figure 90: Manual lamp change

1	Latch	2	Housing doors
3	Cover	4	Push lever
5	Bulb 1	6	Bulb 2

Prerequisite

- ☒ Remove the power plug from the power outlet to safely disconnect the device from the power supply.
- ☒ Watch out for obstacles when changing lamps: Note the position of the microscope!

Action

1. Grasp the cover of the lamp housing by the recessed grips on the left and right side with both hands.
2. Pull the cover off of the device in a horizontal direction.
3. Open the latch of the housing doors.
4. Open the housing doors.
5. Slide or press the push lever towards the second lamp which is not in use as far as it will go.
6. Close the housing doors and latch them again.

7. Reattach the cover to the receptacles on the device.

⇒ The cover is properly attached if the receptacles engage fully.

8. Plug the power cord back into the power outlet.

9. Be sure to replace the defective lamp with a new one [257] following surgery.

10. Dispose of the old lamp according to your local directives and laws.

9.8 Replacing the lamp container

CAUTION!

Risk of injury and burns!

In case of malfunction, the high pressure inside the hot lamp may cause the lamp to burst. Also, the hot surface of the lamp container may also cause burns.

- ▶ Let the lamp container cool off for at least 10 minutes before replacing it.

CAUTION!

Improper handling of the lamp container may result in injuries!

Improper handling of the xenon lamp can lead to damage or injury.

- ▶ The lamp container may be changed only by properly trained persons.


Figure 91: Lamp container, replacement of

1	Latch	2	Housing doors
3	Cover	4	Old lamp container
5	New lamp container		

Prerequisite

- ☑ Remove the power plug from the power outlet to safely disconnect the device from the power supply.
- ☑ Watch out for obstacles when changing lamps: Note the position of the microscope!

Action

1. Grasp the cover of the lamp housing by the recessed grips on the left and right side with both hands.
2. Pull the cover off of the device in a horizontal direction.
3. Open the latch of the housing doors.
4. Open the housing doors.
5. Pull out the defective lamp container.
6. Slide a new lamp container into the lamp housing as far as it will go.
7. Close the housing doors and latch them again.
8. Reattach the cover to the receptacles on the device.
 - ⇒ The cover is properly attached if the receptacles engage fully.
9. Plug the power cord back into the power outlet.

Following a lamp failure, a message stating that you should keep an extra new lamp container on hand as a replacement appears on the monitor during every device start.

Empty page, for your notes

10 Technical specifications

10.1 Regulatory information

The ZEISS KINEVO 900 fulfills the requirements of the following standards:

- IEC 60601-1
- IEC 60601-1-2
- IEC 60825
- IEC 62304

Classification of the device in accordance with IEC 60601

The ZEISS KINEVO 900 is classified as follows:

- Degree of protection against electric shock: Class 1
- Operating mode: continuous operation

Classification of the device in accordance with IEC 60601-1-2

- Electromagnetic compatibility (EMC): Fulfills IEC 60601-1-2, Class A (as per CISPR 11)

10.2 Radio Frequency Approval

This device contains radio modules and fulfills the requirements of 2014/53/EC.

10.3 Radio modules

10.3.1 Bluetooth module

FCP WL, FCP Gateway WL

Designation	Value
Transmit and receive frequencies	2402 MHz to 2480 MHz
Receiving power	-82 dBm to 0 dBm
Transmission power	1 mW to max. 2.5 mW (Class 2)
Modulation	FHSS

10.3.2 WLAN module (option)

The WLAN module cannot be used in all countries.

Intel® 7260HMW

Designation	Value
Transmit and receive frequencies	2412 MHz to 2472 MHz
Transmission power	The maximum output power may vary according to country depending on local regulations.*
Modulation	CCK, DQPSK, DBPSK, BPSK, QPSK, 16 QAM, 64 QAM, 256 QAM
Safety standard	WPA2/PSK

* Output power


- 802.11b: +16 dBm minimum
- 802.11g: +14 dBm minimum
- 802.11a: +14 dBm minimum
- 802.11n HT20 (2.4 GHz): +13 dBm minimum
- 802.11n HT40 (2.4 GHz): +13 dBm minimum

10.3.3 RFID module

TWN4 MultiTech

Designation	Value
Transmit and receive frequencies	13.56 MHz
Transmission power	-54 dBm
Modulation	ASK

10.4 Electrical data

Name	Value
Rated voltage	100 - 240 V AC
Current consumption	max. 1350 VA
Rated Frequency	50 - 60 Hz
Degree of protection (as per EN 60601-1)	I
Fuse	Automatic circuit breaker
USB 3.0 connector	Rated value: 5 V Permissible: 4.45 - 5.5 V / max. 0.9 A*
Remote connector (AUX) 	Max. 24 V / 0.5 A, galvanically decoupled as per IEC 60601-1 Pin assignment: Pin 1: NC Pin 2: NO PIN 3: C

* A device may be supplied with this current intensity only following release by the host controller. Until then, a maximum current intensity of 0.1 A applies.

10.5 Light source

Name	Value
Color temperature	5000 K (±500 K)
Rated power	Approx. 300 W

Additional information

- Illumination technique: Fiber optics
- Main bulb: Xenon short-arc reflector lamp,
- Spare bulb: Xenon short-arc reflector lamp
- Bulb change: Automatic / manual to spare bulb
- Replacement of lamp container: Manual
- Filter: UV/IR heat protection filter, sieve aperture, controls light intensity

10.6 Aiming beam laser

Name	Value
Max. power (IEC 60825-1)	< 1 mW
Wavelength	λ 635 - 645 nm
Laser class	II

Function of laser beam

The focal point is determined by two visible laser beams which intersect in the focal plane.

10.7 Integrated HD camera

Designation	Value
Resolution	1920 x 1080 p
Signal-to-noise ratio	54 dB
Sampling frequency	50 Hz, 59.94 frames/second

10.8 4K camera

Designation	Value
Resolution	3840 x 2160 p
Signal-to-noise ratio	54 dB
Sampling frequency	50 Hz, 59.94 frames/second

10.9 Digital video outputs

Designation	Value
Output: 6 HDMI / DVI	1920x1080p50/60
Output: 7, 8 Display port	1920x1080p50/60
Output: 9, 10 HD-SDI (-0.8 Vp-p/75 Ω PAL)	1920x1080p50/60

10.10 Digital video inputs

Designation	Value
Input 3, 4: Display port	1920x1080p50/60
Input: 5 HDM / DVI	1920x1080p50/60
	1280x720p50/60

10.11 4K video outputs

Name	Value
Outputs 1 and 2 HDMI UHD 4K	HDMI 2.0 3840x2160 50/60p (2xODU)
Color sampling	4:2:0
Outputs 11 to 16 3G-SDI	QuadSDI 3840x2160 50/60p (4xBNC)

10.12 Audio input

Designation	Value
Plug-in connection	3.5 mm stereo jack plug
Voltage level	-10 dBV

Additional information

- When connecting microphones, the signal must be strengthened in advance.
- The microscope-internal microphone must be deactivated when connecting microphones to avoid disruptive noises.

10.13 Mechanical data

10.13.1 System data

Name	Value
Max. additional load on microscope body	Max. 6 kg
Max. dimensions in transport position (WxHxD)	850 x 1900 x1700 mm
Total weight of device with max. add. load	Max. 395 kg
Total weight of device incl. single- use transport box	Approx. 525 kg

10.13.2 Dimensions of rollers

Name	Value
Diameter	150 mm
Tire width	15 mm
Roller width	70 mm

10.13.3 Component weights

The information in this paragraph only applies to the optional configuration of the microscope as a system for hybrid (optical-digital) visualization.

Name	Weight
Stereo co-observation module (000000-1063-869)	1084 g
Tiltable tube, 180° (303791-0000-000)	860 g
Straight tube, f = 170 mm (303765-0000-000)	528 g
Foldable tube, f = 170/260 mm (303771-9021-000)	940 g
Folding tube, f = 170 / 260 mm for mouth switch (303771-9110-000)	940 g
Angle optics with dovetail (spine adapter) (302581-9200-000)	440 g
Photo adapter f= 340 mm (000000-1022-973)	420 g
Photo adapter, T2 Canon EOS (000000-0448-028)	60 g
Mouth switch for 180° tiltable tube (000000-1177-805)	470 g
Push-in widefield eyepieces, 10x (2 pcs) (305542-0000-000)	216 g
Push-in widefield eyepieces, 12.5x (2 pcs) (305543-9901-000)	216 g
Magnification changer, 3-position (303429-9903-000)	448 g
Rotating adapter (301007-0000-000)	192 g

Name	Weight
Dovetail guide (303360-9903-000)	40 g
Micromanipulator	See manufacturer documentation, max. 1 kg.

10.14 Optical data

10.14.1 Surgical microscope

Name	Value
Varioskop	approx. 200 ... 625 mm
Zoom system magnification factor	0.4x - 2.4x
Aiming beam laser (Autofocus option):	
Max. power (IEC 60825-1)	< 1 mW
Shaft length	λ 635 - 645 nm

Additional information

- Magnification adjustment: motorized/manual
- Focusing: motorized/manual
- Sensor for working distance and magnification for neuronavigation
- Motor-driven, travel speed adjustable, automatic adaptation to magnification

10.14.2 Widefield eyepiece (magnification factor 10x)

Designation	Value
Focal Length	25 mm
Field of view	21 mm
Distance of the exit pupil from the last lens	24 - 25.5 mm
Diopter adjustment range	+5/-8
Weight	120 g

10.14.3 Widefield eyepiece (magnification factor 12.5x)

Designation	Value
Focal length	20 mm
Field of view	18 mm
Distance of the exit pupil from the last lens	22 - 23.5 mm
Diopter adjustment range	+5/-8
Weight	115 g

10.15 Ambient requirements for operation

Name	Permissible range
Temperature	+10 ... +40°C
Relative humidity	30 ... 75%
Air pressure	700 ... 1060 hPa
Slope of ground/floor (normal firm ground, solid floor)	Max. 1.5°

10.16 Ambient requirements for transport and storage

Designation	Permissible range
Temperature	-20 ... +60 °C
Rel. humidity (without condensation)	10 ... 90%
Atmospheric pressure:	500 ... 1060 hPa

10.17 Operation with intraoperative MRT systems

For operation in the vicinity of an MRT/MRI scanner, the device must be positioned so that all system components (stand, microscope) are located outside of the 5-Gauss line. The operator must take suitable measures to ensure that the device cannot be moved toward the MRT/MRI system across the 5-Gauss line. The device must be switched off during an MR image acquisition. You may use the device only when no MR image acquisition is in progress.

10.18 Dimensional drawing


Figure 92: Dimensional drawing (dimensions in mm)

10.19 Guidelines and manufacturer's declaration for electromagnetic compatibility

10.19.1 EMC - Electromagnetic compatibility as per IEC 60601-1-2: 2007 (3rd Edition)

The device is subject to specific precautions with regard to electromagnetic compatibility (EMC). In order to avoid the occurrence of EMC interference, the device may only be installed, operated and maintained in the manner indicated in these Instructions for Use and only with components supplied by ZEISS.

NOTE

Danger from electromagnetic radiation!

The ZEISS KINEVO 900 may be disturbed by other devices even when these other devices comply with the emission requirements applicable to them according to CISPR.

- ▶ Do not use the ZEISS KINEVO 900 when it is located next to or stacked on top of other devices.
- ▶ If operation of the device located next to or stacked on top of other devices is required, observe the ZEISS KINEVO 900 to ensure its normal operation in the arrangement in which it is used.

NOTE

Danger from electromagnetic radiation!

Electrical devices can influence each other as a result of their electromagnetic radiation. The use of non-approved components (accessories, transformers of all types, cables) can cause increased emissions or reduce the device's immunity.

- ▶ Only use accessories, transformers, cables and spare parts which are specified in these Instructions for Use or which are approved by ZEISS for this device.
- ▶ Do not use any portable or mobile RF communication equipment near the device as it is not possible to exclude the possibility that the function of the device will be affected.
- ▶ Please follow the EMC guidelines in the following pages.

10.19.1.1 Electromagnetic interference

ZEISS KINEVO 900The is intended for operation in an electromagnetic environment as specified below. The customer or the user of the ZEISS KINEVO 900 is responsible for ensuring that the device is operated in such an environment.		
Interference measurements	Compliance	Electromagnetic environment - guidelines
RF emissions as per CISPR 11	Group 1	The ZEISS KINEVO 900 uses RF energy only for its internal functions. As a result, RF emissions are very low and unlikely to cause any interference in nearby electronic devices.
RF emissions as per CISPR 11	Class A	ZEISS KINEVO 900 is suitable for use in facilities other than residential environments that are directly connected to a PUBLIC POWER GRID which also supplies buildings used for residential purposes.
Harmonic emissions as per IEC 61000-3-2	Not applicable	
Emission of voltage fluctuations/ flicker as per IEC 61000-3-3	Not applicable	
NOTE The properties of the ZEISS KINEVO 900 determined by EMISSIONS permit its use in industrial areas and in hospitals (CISPR 11, Class A). If used in a residential area (for which Class B usually is required as per CISPR 11), the ZEISS KINEVO 900 may not provide the necessary protection with regard to radio frequency communication services. The user may have to take corrective measures such as moving or realigning the device.		

10.19.1.2 Electromagnetic immunity for all ME equipment and ME systems

ZEISS KINEVO 900The is intended for operation in an electromagnetic environment as specified below. The customer or the user of the ZEISS KINEVO 900 is responsible for ensuring that the device is operated in such an environment.			
Electromagnetic immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment - guidelines
Electrostatic discharge (ESD) as per IEC 61000-4-2	±6 kV contact discharge ±8 kV air discharge	±6 kV contact discharge ±8 kV air discharge	Floors should be made of wood or concrete or be covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Fast transient/ burst immunity as per IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	The quality of the supply voltage should be that of a typical business or hospital environment.

Surges as per IEC 61000-4-5	± 1 kV line-to-line voltage ± 2 kV line-to-ground voltage	± 1 kV line-to-line voltage ± 2 kV line-to-ground voltage	The quality of the supply voltage should be that of a typical business or hospital environment.
Voltage dips, short interruptions and voltage variations as per IEC 61000-4-11	$< 5\% U_T$ ($> 95\%$ dip in U_T) for 1/2 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $< 5\% U_T$ (95% dip in U_T) for 5 s	$< 5\% U_T$ ($> 95\%$ dip in U_T) for 1/2 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $< 5\% U_T$ (95% dip in U_T) for 5 s	The quality of the supply voltage should be that of a typical business or hospital environment. If the user of the ZEISS KINEVO 900 requires continued operation even in the event of interruptions in the power supply, we recommend powering the ZEISS KINEVO 900 from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field as per IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields in the supply frequency should correspond to the typical values that are found in business and hospital environments.

Note: U_T is the AC voltage supply before application of the test levels.

10.19.1.3 Electromagnetic immunity for non-life-supporting devices

The ZEISS KINEVO 900 is intended for operation in the electromagnetic environment specified below. The customer or the user of ZEISS KINEVO 900 is responsible for ensuring that the device is operated in such an environment.			
Electromagnetic immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment - guidelines
Conducted RF disturbances as per IEC 61000-4-6 Emission of RF disturbances according to IEC 61000-4-3	$3 \text{ V}_{\text{effective value}}$ 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 V/m 3 V	Portable and mobile radio communication equipment should not be used closer to ZEISS KINEVO 900, including its cables, than the recommended safety distance that is calculated using the equation applicable to the transmission frequency involved. Recommended safety distance $d = 1.2 \sqrt{P}$

			$d = 1.2 \sqrt{P}$ for 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ for 800 MHz to 2.5 GHz where P is the rated output of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). The field strength of stationary transmitters as determined by a site survey^a should be less than the compliance level in all frequency ranges.^b Interference may occur in the vicinity of equipment marked with the following symbol. 
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is influenced by absorption and reflection by structures, objects and persons.			
^a Theoretically, field strengths of stationary transmitters such as base stations for mobile telephones and mobile land radio equipment, amateur radio stations, AM and FM radio broadcast and TV broadcast transmitters cannot be predicted accurately. To assess the electromagnetic environment with respect to stationary RF transmitters, a site study of the electromagnetic phenomena should be considered. If the measured field strength in the location where the ZEISS KINEVO 900 is used exceeds the COMPLIANCE LEVELS indicated above, the ZEISS KINEVO 900 should be monitored to verify normal OPERATION. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the ZEISS KINEVO 900. ^b Field strengths should be less than 3 V/m over the frequency range from 150 kHz to 80 MHz.			

10.19.1.4 Recommended protective distances between bearable and mobile RF telecommunications systems and the ME equipment

ZEISS KINEVO 900 is intended for use in an ELECTROMAGNETIC ENVIRONMENT in which RF disturbances are controlled. The customer or the user of the ZEISS KINEVO 900 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the ZEISS KINEVO 900 — depending on the output power of the communication equipment as specified below.

Rated output power of the transmitter W	Safety distance, dependent on transmission frequency m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1.2 \sqrt{P}$	$d = 1.2 \sqrt{P}$	$d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation indicated for each column, with P being the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer's specifications.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is influenced by absorption and reflection by structures, objects and persons.

10.19.2 EMC - Electromagnetic compatibility as per IEC 60601-1-2: 2014 (4th Edition)

The device is subject to specific precautions with regard to electromagnetic compatibility (EMC) in Professional Healthcare Facility Environments.

In order to avoid the occurrence of EMC interference, the device may only be installed, operated and maintained in the manner indicated in these Instructions for Use and only with components supplied by ZEISS.

WARNING!

Function deterioration!

Do not install or operate the device in direct proximity to other devices with the exception of the combination of the devices described in these Instructions for Use.

- If it cannot be avoided that the ZEISS KINEVO 900 is operated in proximity to other devices, the proper function of the ZEISS KINEVO 900 must be monitored.

WARNING!

Function deterioration!

- ▶ Only use accessories, transformers, cables and spare parts which are specified in these Instructions for Use or which are approved by ZEISS for this device.

WARNING!

Deterioration of performance!

- ▶ Do not use any portable or mobile HF communication equipment or transmitters (including peripheral devices such as antenna cables or external antennas) in the proximity of the device (minimum distance 30 cm), as it cannot be ruled out that the function of the device will be impaired or that the performance of the device deteriorates.
- ▶ Please follow the EMC guidelines in the following pages.

NOTE

Danger from electromagnetic radiation!

The ZEISS KINEVO 900 may be disturbed by other devices even when these other devices comply with the emission requirements applicable to them according to CISPR.

- ▶ Do not use the ZEISS KINEVO 900 when it is located next to or stacked on top of other devices.
- ▶ If operation of the device located next to or stacked on top of other devices is required, observe the ZEISS KINEVO 900 to ensure its normal operation in the arrangement in which it is used.

10.19.2.1 Electromagnetic interference

ZEISS KINEVO 900 is intended for operation in an electromagnetic environment as specified below. The customer or the user of the ZEISS KINEVO 900 is responsible for ensuring that the device is operated in such an environment.

Interference measurements	Compliance
RF emissions as per CISPR 11	Group 1
RF emissions as per CISPR 11	Class A
Emission of voltage fluctuations/flicker as per IEC 61000-3-3	Not applicable

NOTE

The properties of the ZEISS KINEVO 900 determined by EMISSIONS permit its use in industrial areas and in hospitals (CISPR 11, Class A). If used in a residential area (for which Class B usually is required as per CISPR 11), the ZEISS KINEVO 900 may not provide the necessary protection with regard to radio frequency communication services. The user may have to take corrective measures such as moving or realigning the device.

NOTE

The emission properties of the ZEISS KINEVO 900 are suitable for use in industrial areas and hospitals (CISPR 11 Class A). If the ZEISS KINEVO 900 is used in a residential area (for which CISPR 11 Class B is normally required) this equipment may not provide the necessary protection with regard to radio frequency communication services.

10.19.2.2 Electromagnetic immunity for all ME equipment and ME systems

ZEISS KINEVO 900The is intended for operation in an electromagnetic environment as specified below. The customer or the user of the ZEISS KINEVO 900 is responsible for ensuring that the device is operated in such an environment.		
Electromagnetic immunity tests	IEC 60601 test level	Compliance level
Electrostatic dis- charge (ESD) as per IEC 61000-4-2	±8 kV contact discharge ±15 kV air discharge	±8 kV contact discharge ±15 kV air discharge
Fast transient/ burst immunity as per IEC 61000-4-4	±2 kV for power lines ±1 kV for input/output lines	±2 kV for power lines ±1 kV for input/output lines
Surges as per IEC 61000-4-5	±1 kV voltage Phase neutral conductor ±2 kV voltage Phase neutral conductor-earth	±1 kV voltage Phase neutral conductor ±2 kV voltage Phase neutral conductor-earth
Power frequency (50/60Hz) mag- netic field as per IEC 61000-4-8	30 A/m	30 A/m
Voltage dips, short inter- ruptions and voltage variations as per IEC 61000-4-11	0 % U_T for 1/2 cycle	0 % U_T for 1/2 cycle
	0 % U_T for 1 cycle	0 % U_T for 1 cycle
	70 % U_T for 25/30 cycles	70 % U_T for 25/30 cycles
	0 % U_T for 250/300 cycles	0 % U_T for 250/300 cycles

10.19.2.3 Electromagnetic immunity for non-life-supporting devices

The ZEISS KINEVO 900 is intended for operation in the electromagnetic environment specified below. The customer or the user of the ZEISS KINEVO 900 is responsible for ensuring that the device is operated in such an environment.		
Electromagnetic immunity tests	IEC 60601 test level	Compliance level
Conducted RF disturbances as per IEC 61000-4-6	3 V 150 kHz to 80 MHz	3 V
	6 V ISM bands between 150 kHz and 80 MHz	6 V
Emission of RF disturbances according to IEC 61000-4-3	3 V/m 80 MHz – 2.7 GHz	3 V/m
Radiated RF disturbances from near fields of wireless communication devices as per EN 61000-4-3	27 V/m 385 MHz	27 V/m
	28 V/m 450 MHz, 810 MHz – 2.45 GHz	28 V/m
	9 V/m 710 MHz – 780 MHz, 5.24 GHz – 5.785 GHz	9 V/m

10.20 Security Concept

10.20.1 Objective

The objective of this document is to summarize the set of cybersecurity controls to ensure medical device cybersecurity for the ZEISS KINEVO 900 of ZEISS with an embedded Windows 10 IoT Enterprise operating system.

10.20.2 System overview

The device has the following external interfaces, which are critical for cybersecurity:

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10.20.3 General principles

Medical device security is a shared responsibility between stakeholders, which include healthcare facilities, patients, service providers and manufacturers of medical devices.

Failure to maintain cybersecurity can result in compromised device functionality, loss of data (medical or personal) availability or integrity, or expose other connected devices or networks to security threats.

10.20.3.1 Identification of general security threats and vulnerabilities

- Malware could repurpose the system's capabilities if it accesses the underlying operating system.
- Malware could change system capabilities and even render the whole system unusable.
- Malware could expose the device to unauthorized use or alter it.
- Unauthorized access could expose the device to unauthorized use or alter it.
- Unauthorized access could compromise data.
- Operators could be tempted to repurpose the system's computing capabilities.

10.20.3.2 Medical device regulations and guidelines

- The medical device as well as any subsequent design changes will be validated and documented by considering the actual regulations and guidelines for medical devices.
- To ensure the proper functionality of the device and to avoid any non-conformities regarding medical device regulations, neither third-party software (e.g. third-party virus-scanning software) nor automatic updates (e.g. automatic Windows updates) may be applied to the system. The same applies to non-validated operating system configuration changes.
- An assessment regarding the impact of cybersecurity threats and vulnerabilities regarding device functionality as well as end users/patients will be performed.
- If threats and vulnerabilities are exploited, the likelihood and risk levels will be determined.
- Based on the respective assessments, subsequent actions to maintain the intended system functionality and security will be taken by ZEISS as and when needed.

10.20.4 Cybersecurity relevant functions

10.20.4.1 Limit access to trusted users only

10.20.4.1.1 User authentication

- User authentication can be enabled (optional).
- For privileged users, password authentication is required.
- Default user for emergency cases is always accessible with restricted functionality (access only to active patient).

10.20.4.1.2 Layered authorization model based on the user role

A layered authorization model with different privileges based on the user role is implemented:

- Standard users have standard basic rights.
- The "IT Admin" has additional IT system rights to enable/set passwords for standard users, to enable/set up network connections (such as DICOM, Remote Service, WLAN, Hotspot and WebUI) and to set the time and date.
- In addition to the "IT Admin", "Biomed" has limited service functionalities available.
- The "ZEISS Service" user has additional rights to access service dialogs and settings.
- The "ZEISS Factory" user has full access.

10.20.4.1.3 Appropriate authentication

With activated user authentication, authentication is required for all users (except default users in cases of emergency):

- Password must be entered on a virtual keyboard displayed on the GUI.
- In the ZEISS KINEVO 900 software, no user authentication is possible via the network.
- Brute-force attacks require time for entering passwords through virtual keyboard and GUI.

10.20.4.1.4 Strong password protection

"IT Admin", "Biomed", "ZEISS Service" and "ZEISS Factory" user authentication passwords are more than 12 characters long.

10.20.4.1.5 BIOS

BIOS is password-protected.

10.20.4.1.6 User authentication for software or firmware updates

Software and firmware updates require the BIOS password.

10.20.4.2 External interfaces

Only interfaces which are necessary for the intended use are accessible. All other standard PC interfaces are located inside the device.

10.20.4.2.1 USB drive file system

- Only data exchange via external storage media
- Controlled via ZEISS KINEVO 900 device software
- Only patient data export/import
- Settings export/import
- License import
- Autorun functionality disabled
- Keyboards via USB not active by default

10.20.4.2.2 RS232 Interface

- Proprietary interface for navigation interface
- No access to patient data
- Just controlled by application
- No access to device application

10.20.4.2.3 Network-based functionalities and network interface bindings

The purpose of the medical device's connection to a computer/IT network:

- Patient data (images, video) export over network to central archiving system via DICOM (LAN, WLAN)
- Patient data (images, video) export over network to shared folder for archiving data (LAN, WLAN)
- Patient data (images, video) import over network from shared folder (LAN, WLAN)
- Streaming of camera's live feed over network to a compatible video streaming client (LAN, WLAN, WLAN hotspot)
- Parallel video recording over network to folder share (LAN, WLAN)
- ZEISS Smart Services integrated remote support client - see also the White Paper ZEISS Smart Services (LAN, WLAN)
- Web browser access from mobile device to WebUI (WLAN hotspot)

- Interface for IOS apps (WLAN hotspot)
- Navigation interface (LAN)

10.20.4.2.4 Network specifications for the respective network interface LAN

- TCP/IP IPv4 Ethernet Protocol
- RTP for video streaming services
- Gigabit-Ethernet
- Static IP or DHCP (configurable within the ZEISS KINEVO 900 system)
- DICOM protocol (no encryption)
- Remote network share access via SMB
- Navigation: SSL encrypted
- ZEISS Smart Services: SSL encrypted

WLAN

- WPA2/PSK
- WPA2/PEAP (WLAN Enterprise)
- 802.11a/b/g/n
- TCP/IP IPv4 Ethernet Protocol
- DICOM protocol (no encryption)
- RTP for video streaming services
- Static IP or DHCP (configurable within the ZEISS KINEVO 900 system)
- Remote network share access via SMB
- ZEISS Smart Services: SSL encrypted

WLAN-Hotspot

- WPA2/PSK
- 802.11a/b/g/n
- TCP/IP IPv4 Ethernet Protocol
- Hotspot by default deactivated, manual activation for use required
- RTP for video streaming services
- DHCP server
- REST API for IOS apps via HTTP (authentication mandatory)
- WebUI via HTTP (authentication mandatory, not encrypted)

10.20.4.3 Firewall

- The system is secured via a built-in Windows firewall.
- The ZEISS KINEVO 900 software is allowed to bind arbitrary ports in- and outbound. However, the ZEISS KINEVO 900 software cannot change the firewall configuration.

The following applications and services are activated and enabled through the built-in firewall:

- DHCP client
- DHCP server (only for WLAN hotspot)
- DICOM service with configurable port
- SMB protocol TCP port 445 (outbound only)
- Navigation interface: TCP port 8123
- RTP video streaming: TCP port 8554
- Ping
- Tracert
- WebUI TCP Ports 80 (HTTP) and 8443 (HTTPS)
- REST API for apps: TCP ports 80 (HTTP) and 8443 (HTTPS)
- ZEISS Smart Services:
 - Remote support functionality
 - TCP outbound, remote destination port 7778
 - Device lookup, inbound ports: UDP 7716–7719
 - Device lookup outbound ports: UDP remote ports 7700 – 7703

The lookup ports will only be used during the automatic, initial setup of the connection between the ZEISS KINEVO 900 and the remote support site control unit (SCU). Alternatively, manual configuration of the connection is possible, so the UDP protocol is not required for remote support functionalities.

All unnecessary services and applications are deactivated by default. Access to the services is controlled by the ZEISS KINEVO 900 software and will only be activated on request

10.20.4.4 Features that allow security compromises to be detected and qualified

- Audit trail, security logs, error logs and user logs are implemented.
- Log files can be exported (encrypted) and sent to ZEISS.
- No personal data (user data or patient data) in log files.
- A on-time virus scan can be performed by a ZEISS service technician in justified, exceptional cases.

10.20.4.5 Protection of Windows operating system inside a medical device

- Customized Windows 10 IoT Enterprise operating system. Unnecessary components are removed.
- Closed system (no direct access to the Windows operating system):
 - No possibility for either users or IT administrators to gain access to the Windows operating system GUI, on which the ZEISS KINEVO 900 software is based
 - Therefore, there is no possibility to install additional software and/or change or modify existing software modules which might introduce any viruses or malware to the system or damage it.
 - It is not possible to execute any files on the system which might act as malware as the user and IT administrator have no access to the Windows operating system GUI.
 - The ZEISS KINEVO 900 software only grants access through the ZEISS KINEVO 900 GUI of the system.
 - Application in full-screen mode. No Windows desktop
 - No hotkeys are accepted by the system as any keyboards are blocked. In addition, system-critical tools like Windows Task Manager are blocked by policies on the system.
- Hard disk (SSD) partitions, where the operating system and ZEISS KINEVO 900 software is stored, are protected against persistent changes (any unauthorized software couldn't be installed persistently to the SSD partition and would disappear after a system reboot; the same applies for files and configuration changes which have been applied to the system)
- BIOS is password-protected.
- BIOS settings only allow booting from the dedicated internal hard disk (SSD). Booting from external devices is disabled. Enabling other boot devices requires a BIOS password.
- The system features a keyboard lock that prevents unauthorized users from launching a Windows Explorer window from a connected USB keyboard. Furthermore it is not possible to execute any third-party software application on a connected external USB device which might act as malware, as only ZEISS field service engineers have access to the Windows 10 operating system.
- The ZEISS KINEVO 900 features a specific HDD design:
 - The Windows 10 operating system and the ZEISS KINEVO 900 software are located on an SSD partition with protection against persistent changes to prevent the installation of unauthorized third-party software applications. Any potential changes by a virus to that partition will be discarded automatically when the system is restarted.

- All user data are stored on the second SSD partition with no write protection.
- Patient data are located on a separate HDD that can be replaced in case of a hardware malfunction.
- Log off and switch user functions are disabled.

10.20.5 Ensure trusted system setup

- Restrict software or firmware updates to authenticated code.
- Software and firmware updates may only be performed by ZEISS field service engineers with specially configured ZEISS software update media.
- All updates require the BIOS password.
- Software checks firmware versions for compatibility.
- Use systematic procedures for authorized users to download version-identifiable software and firmware from the manufacturer.
- No software downloads are provided to customers (only to ZEISS field service engineers).
- During manufacturing and after each configuration change by a ZEISS field service engineer, a backup is performed and stored.
- ZEISS field service engineers can restore any stored backup.
- Images for the ZEISS KINEVO 900 software and software updates are checked for viruses.
- To ensure a verified medical device setup and functionality and to maintain medical device conformity, no non-validated changes to the system are allowed:
 - The automatic Windows update function is deactivated by design for the above reason. However, the ZEISS KINEVO 900 system in its hardened configuration offers no target for the exploitation of any system vulnerabilities. None of the ZEISS KINEVO 900 or system services are represented by default in an uncontrolled manner to the outside of the system.
 - For the same reason, no third-party virus-scanning solutions may be installed and run on the ZEISS KINEVO 900 system. To protect the system against viruses and other malware, all options to execute any files or applications are controlled by the ZEISS KINEVO 900 software. There is no way to execute files or applications either manually or automatically on the ZEISS KINEVO 900 system. As any changes to the partitions on which the operating system and the ZEISS KINEVO 900 software reside are non-persistent, there is no way for any malware to be installed and executed on a ZEISS KINEVO 900 system.

- The same applies to linking a ZEISS KINEVO 900 to an AD domain.
- Windows operating system patches will be delivered and installed by ZEISS only, as and when required.

10.20.6 Protection of critical device functionality

The device magnification, focus adjustment, light intensity and manual ZEISS KINEVO 900 movement features are independent of any PC functionality. The mentioned device features can still be used even if the integrated PC is not available or shut down (basic function mode).

Empty page, for your notes

11 Accessories and components

These instructions for use describe accessories that are not essential components of the individual deliveries. A current accessories and components list can be obtained from your ZEISS contact person.

You can find the ZEISS contact partner for your country on the following website: www.zeiss.com/med

Only use accessories and components authorized by ZEISS for this device. Device safety during operation cannot be guaranteed if accessories and components that are not authorized by ZEISS are used.

11.1 Video components

Video components can be found in the separate product overview G-30-1888.

11.2 QEVO

Designation	Specification	Order no.
QEVO	Digital Exploration Tool	303155-9020-000

11.3 Components for hybrid visualization

The components specified in this paragraph are only applicable if the microscope is configured as a system for hybrid (optical-digital) visualization.

11.3.1 Main tube

Designation	Specification	Order no.
Tiltable tube 180°	Swivel range 180°, f=170mm	303791-0000-000
Foldable tube, white	f=170/260 mm not suitable for mouth switch adaptation	303771-9021-000
Foldable tube for mouth switch, white	f=170/260 mm without integrated rotation function (only in conjunction with mouth switch 000000-1177-805)	303771-9111-000
Angle optics with dovetail (spine adapter), white		302584-9200-000

11.3.2 Eyepieces for main tube

Designation	Specification	Order no.
Push-in widefield eyepieces, 2 pcs	10x	305542-0000-000
Push-in widefield eyepieces, 2 pcs	12.5x aspherical	305543-9901-000

11.3.3 Co-observer tube

Designation	Specification	Order no.
Stereo co-observation module	with 2 pivot joints and locking levers	000000-1063-869

11.3.4 Tube for left/right co-observation

Designation	Specification	Order no.
Straight tube	f=170 mm	303765-0000-000
Sleeves for tubes with screw thread	2x	305542-0107-000
180° tiltable tube	Swivel range 180°, f=170mm	303791-0000-000
Foldable tube, white	f=170/260 mm not suitable for mouth switch adaptation	303771-9021-000

11.3.5 Eyepieces for left/right co-observation

Designation	Specification	Order no.
Push-in widefield eyepieces, 2 pcs	10x	305542-0000-000
Push-in widefield eyepieces, 2 pcs	12.5x aspherical	305543-9901-000

11.3.6 Tube for face-to-face

Designation	Specification	Order no.
180° tiltable tube	Swivel range 180°, f=170mm	303791-0000-000
Foldable tube, white	f=170/260 mm not suitable for mouth switch adaptation	303771-9021-000
Rotary adapter	-	301007-0000-000

11.3.7 Eyepieces for face-to-face

Designation	Specification	Order no.
Push-in widefield eyepieces, 2 pcs	10x	305542-0000-000
Push-in widefield eyepieces, 2 pcs	12.5x aspherical	305543-9901-000

11.3.8 Tube for posterior fossa

Designation	Specification	Order no.
Straight tube	f=170 mm	303765-0000-000
Sleeves for tubes with screw thread	2x	305542-0107-000

11.3.9 Eyepieces for posterior fossa

Designation	Specification	Order no.
Push-in widefield eyepieces, 2 pcs	10x	305542-0000-000
Push-in widefield eyepieces, 2 pcs	12.5x aspherical	305543-9901-000

11.3.10 Canon EOS DSLR Adapter Kit

Designation	Specification	Order no.
Photo adaptor	f=340	000000-1022-973
T2 adapter	Canon EOS 30743	000000-0448-028
Control cable	Canon EOS D60	000000-1229-877

11.3.11 Mouth switch

Designation	Specification	Order no.
Mouth switch	for 180° tiltable tube and foldable tube (303771-9110-000)	000000-1177-805
Angle adapter for mouth switch		302700-8602-000

11.3.12 Magnification Changer

Designation	Specification	Order no.
Magnification Changer	3-position	303429-9903-000

11.3.13 Adaptation of laser micromanipulators

Designation	Specification	Order no.
Dovetail mount / adapter plate for ext. components	(e.g. laser micromanipulator)	303360-9903-000

11.4 Tube covers

Designation	Specification	Order no.
Cover, main observer tube	NCS 1002-B	302584-1040-000
Cover, co-observer tube	NCS 1002-B, lak	302584-1041-000

11.5 Fluorescence targets

Designation	Specification	Order no.
Fluorescence target	BL 400	302581-9052-000
Fluorescence target	YE 560	302582-9208-000

11.6 Foot control panel

Designation	Specification	Order no.
FCP WL	Foot control panel with 14 functions, wireless	304970-9200-000
FCP	Foot control panel with 14 functions, cable-connected	304970-9100-000
Foot switch components		
Battery set	3 pieces, packaged	304970-8821-000
FCP indicator label	FCP indicator label	304970-1010-000
Instructions for Use for FCP and FCP WL	G-30-1706 FCP and FCP WL	000000-1520-536

11.7 FCP cable/FCP backup cable

Designation	Specification	Order no.
FCP CAN bus cable	6m	304970-8760-000

11.8 Rocker foot switch

Designation	Specification	Order no.
Rocker switch	2 functions, 3 m cable	305989-8609-000

11.9 Sterile covers/drapes

Designation	Specification	Order no.
SMARTDRAPE	No. 28 (5 pcs), sterile	306028-0000-000
VisionGuard replacement lenses	(20/box), sterile	306001-0000-000

11.10 Country-specific cables

Designation	Specification	Order no.
Europe	Power cord, length: 6 m	000000-0594-821
Netherlands	Power cord, length: 6 m	000000-0603-410
USA	Power cord, length: 6 m	000000-0594-822
UK	Power cord, length: 6 m	000000-0594-823
Switzerland	Power cord, length: 6 m	000000-0584-947
Argentina	Power cord, length: 6 m	000000-0594-906
China	Power cord, length: 6 m	000000-0594-824
Brazil	Power cord, length: 6 m	000000-0594-905
Australia / New Zealand	Power cord, length: 6 m	000000-0616-997

11.11 Stereo glasses

Designation	Specification	Order no.
Set of stereo polarization glasses	5 pairs	000000-1992-943

11.12 Medical loupes

Designation	Specification	Order no.
EyeMag Pro S	System carrier S (support) with components	304160-9100-000
	2x touch guard (EyeMag Pro F/S)	304156-9900-000
	Large softcase for Pro S	000000-1460-441
	Outer packaging for Pro S	000000-0488-101
	Optical system K 4.3x/500 silver	304156-9404-000

11.13 Lamp container of the xenon light source

Designation	Specification	Order no.
Lamp container	HLQ300W+	304949-9002-000

11.14 USB stick

Designation	Specification	Order no.
USB stick	USB 3.0 Memory Stick 64 GB	000000-0576-392

11.15 USB memory medium

Designation	Specification	Order no.
USB memory medium	External Mini HDD USB 3.0 1 TB	000000-0574-996

12 Disposal

NOTE

Incorrect disposal contaminates the environment!

A high proportion of electronic scrap at the unsorted municipal waste facility poses a risk to the environment.

- ▶ Do not dispose of the waste as normal municipal waste.
- ▶ Observe the local laws/regulations governing the disposal of electrical and electronic equipment.

- ▶ Keep packing material in the event of a relocation or repair.
- ▶ If you would like to dispose of the original packing material, send it in for recycling via a recognized collection system.

The device contains electronic components with integrated batteries.

- ▶ Dispose of the device and integrated batteries correctly, in accordance with national legislation.

In accordance with applicable EU guidelines and national regulations at the time at which the product was brought onto the market, the product specified on the consignment note is not to be disposed of via the domestic waste disposal system or communal waste disposal facilities.

- ▶ For more information on disposing of the device, contact the ZEISS contact person in your country.

You can find the ZEISS contact partner for your country on the following website: www.zeiss.com/med

- ▶ If you resell the device or its components: Inform the buyer that the device is to be disposed of in accordance with the currently applicable regulations.

Empty page, for your notes

Glossary

DHCP

DHCP (Dynamic Host Configuration Protocol) is a network configuration for dynamically assigning IP addresses.

Diopter scale

Element of an eyepiece for reading off the set refraction value.

DNS

Domain name system

Electromagnetic compatibility (EMC)

EMC (electromagnetic compatibility) designates the normal, desired state in which technical devices do not impede each other as the result of undesired electric or electromagnetic effects (non-interference).

Eyecup

An eyepiece control element used to shield the eyepiece against scattered light during eye-controlled focusing.

FL button

Release button for fluorescence application on the handgrip or FCP

IP

The IP (Internet Protocol) is a communications protocol to transfer data in the Internet or a network.

LAN

LAN (local area network)

Ref. Phys. Name

Referring Physician's Name, name of the referring physician who requested the procedure.

Req. Proc. Code

Requested Procedure Code, code value that describes the requested procedure according to a specific coding scheme.

Req. Proc. Desc.

Requested Procedure Description, administrative description specified by the hospital or classification of the requested procedure.

Req. Proc. ID

Requested Procedure ID, identification number that identifies the requested procedure in the "Imaging Service Request".

Sch. Proc. Step Desc.

Scheduled Procedure Step Description, description specified by the hospital or classification of the scheduled procedure step.

Sch. Proc. Step Start

Scheduled Procedure Step Start, planned start time for the procedure step.

Sch. Protocol Code

Scheduled Protocol Code, code value that describes the scheduled protocol according to a specific coding scheme.

SSID

Service Set Identifier, network name of a functional network (e.g. of a WLAN)

Technical safety test

Technical safety test to determine and assess the system safety.

UDI

Unique Device Identification (UDI)
Standardized identification system for medical devices.

UDI Production Identifier (UDI-PI)

Unique Device Identification - Production Identifier

UDI-DI

Unique Device Identification - Device Identifier

USB

USB (universal serial bus) is a standard connector to connect peripheral devices.

WLAN

Wireless Local Area Network

Index

A

AB brake release button	62
AB button.....	181
Acknowledging error messages	58
Activate User Password.....	154
Activate WLAN	151
Adding a patient.....	205
Adding a user	201
Additional Light	55
Additional light on/off	175
Adjust focus/working distance manually.....	60
Adjust luminous field diameter manually.....	59
Adjust pivoting mirror.....	59
Adjust zoom magnification manually ..	60
Adjusting light intensity	175
Anti-fogging agents.....	245
Approval label	32
Attaching a SMARTDRAPE	164
Attaching drapes	164
Attaching tubes	116
Audio configuration menu	88
Auto Brightness	78
Autobalance	165
Automatic circuit breaker.....	46

AUX

Pin assignment.....	263
Remote socket	45

B

Basic Function Mode	254
Binocular tubes	60
BLUE 400 (Option)	232

Button function	64
-----------------------	----

C

Camera release socket	59
Change Date Format.....	155
Change of user	203
Change Time Format.....	155
Changing patients.....	205
Charging status of battery in FCP	70
Cleaning	
Mechanical surfaces	245
Monitor (touchscreen)	244
Optical surfaces.....	244
Clutch release keys.....	181
Component weights	267
Computer fails	254
Configurable buttons	62, 82, 84
Configuration menu displays.....	87
Configure fluorescence button.....	233
Configure network.....	93, 149
Configuring the hand grips	180
Configuring the illumination.....	175
Configuring the Service PC.....	154
Configuring the XY Movement Mode handgrip button	189
Connector panel	45
Control elements for widefield eyepiece	61
Controls of binocular tubes.....	60
Controls on the microscope	59
Coobservation equipment, mounting of.....	120
Co-observation port.....	50
Co-observer (start value)	79

Cranial procedures.....	161	Extras configuration menu	109
Creating a recording.....	217	Eyecup.....	53
Creating photos.....	217	Eyepieces, attaching	116
Current consumption.....	263		
D		F	
Degree of protection	261	Face-to-face (spine) procedures.....	161
Delete User.....	204	Failure of the computer.....	254
Deleting patients	210	FCP (foot control panel), configuration.	181
Device in the de-energized state	255	FCP configuration menu.....	84
Device ON/OFF switch	125	FCP connection status.....	70
DICOM (option)	222	FCP signal strength	70
DICOM network connection	153	FL button	235, 238
Digiskop	118	Fluorescence configuration menu	103
Digiskop cover.....	54, 118	Focal Length.....	79
Dimensional drawing.....	269	Focus configuration menu.....	80
Diopter scale.....	53	Focus Light Link	78
Disinfection	245	Focus Zoom Link.....	81
Disposal		Foot control panel	
Batteries.....	293	Buttons	63
Electronics.....	293	Joystick	63
Packing material.....	293	Rocker switches.....	63
Documentation equipment, mounting of	120	Forgot your password?	149
Drape position	162, 191	G	
E		Green dot	64
Electrical safety	25	H	
EMC malfunctions	21	Hand grip configuration menu	82
Environmental conditions	26	Hand grips, adjustment.....	160
Error message		Help display	57
Export Log Files.....	251	Hotspot	
Export Log Files	251	Connection (display)	70
Export Patient Data.....	207	Display password.....	71
Exporting patient data to DICOM.....	223	Using the device as a hotspot.....	239
Exporting user settings to a USB drive.	202	I	
		Image Rotation	228

Import WLAN certificate	151	Multilayer T* super anti reflection coating	244
Importing patient data.....	208	MultiVision configuration menu	106
Importing patient data via DICOM	209, 223		
Importing user settings from a USB drive	203	N	
Info configuration menu	108	Network drive (network storage)	152
Initial device login.....	127	Non-sterile area	164
Initial instruction.....	15	Normal use	16
Intended use / Indication for use (US only)	15		
Internal hard drive space.....	70	O	
IT admin password forgotten.....	149	Online support.....	57
IT admin system rights	149, 241	Operating mode	261
		Operating personnel	21
		Other applicable documents	13
L			
Labeling on stand	33	P	
Labeling on the microscope	31	Pair wireless foot control panel	129
Lamp change, manual	256	Pairing the FCP WL with the device.....	129
Laser micromanipulators		Parking position	122, 123, 190
Fitting the dovetail mount	157	Password protection for handgrips and FCP	183
Mounting	157	Patient administration	205
Liability	26	Patient video	
Light configuration menu	78	Cutting.....	221
Light Warning Threshold.....	78	Editing.....	221
Light warning threshold adjustment ...	177	Viewing.....	220
Locking mechanism tab	122	Pedal for straight-ahead movement	122
		Photo configuration menu	89
M		PointLock movement mode, manual	186
Magnetic coupling, eyepieces	53	PointLock movement mode, motorized	188
Manage User	201	PointLock, manual	186
Manual lamp change	256	PointLock, motorized	188
Max. additional load on microscope body	265	Potential equalization.....	125
Microscope movement mode	187	Power input socket	125
Mounting the mouth switch	121	Product modification	22
Mouth switch socket	59	Protect Settings	183
Movement in the de-energized state...	255		

Q

QEVO	69, 167,
	168,
	224,
	225, 226
	, 227,
	228, 231

R

Radio frequency approvals	42
Rated frequency	263
Rated voltage	263
Recording configuration menu	102
Recording video files	
In HD quality	210
In LowRes quality	210
Parallel on network drive	213
Parallel on USB	213
Via SmartRecording	214
Remote connector (AUX)	263
Replacing the lamp container	258
RFID reader	164
Risk of burn injuries	18
Rocker foot switch configuration menu	85
Rollers, dimension	266

S

Safety inspection	248
SB brake release button	62
SB Button	181
Searching for a patient	206
Set System Language	155
Setting the date and time	155
Setting the focus and zoom start values	64
Setting up a network connection	
Via DICOM	153

Via LAN	150
Via WLAN	151
SmartLinks configuration menu	77
SmartRecording	214
Sorting patients	206
Stand configuration menu	86
Stand movement mode	188
Standard Password	154
Standby/ON-OFF switch	125, 126
Status Information	70
Sterile area	164
Streaming	216
Switching a device	126
Switch-off	242
System Mode	79, 118

T

Taking a photo with the DSLR camera (option)	218
Target group	15
Transport direction	122
Transport handles	122
Tube configuration menu	79

U

UDI label	32
USB medium storage space	70
USB port	46
User-specific settings	64

V

Video camera	
Menu configuration	194
Menu description	90
Video connector panel	47
video files	
Streaming	216

Video Frequency	155
Video recording	
Cutting.....	221
Editing	221
Parallel Recording.....	213
Viewing.....	220
Viewing images	219
Viewing patient images	219
W	
Warranty	26
Web Interface	
Activation (in the device)	152
Connecting an external device to the KINEVO 900.....	240
Display of address	71
Weight with max. additional load	265
White Balance	196
Wide-field eyepiece	53
X	
XY adjustment mode	184
XY Zoom Link.....	81
Y	
YELLOW 560 (option)	235
Z	
Zoom configuration menu	81

Carl Zeiss Meditec AG

Goeschwitzer Strasse 51-52

07745 Jena

Germany

www.meditec.zeiss.com

Subject to change.

G-30-2094-en - 1.1 - 2020-09-14

FOR THE RUSSIAN FEDERATION ONLY

The Instruction for use «KINEVO 900 Surgical Microscope with accessories», manufactured by "Carl Zeiss Meditec AG" (Carl Zeiss Meditec AG) in accordance with national and international requirements, is standard for all countries. In accordance with national requirements and in order to comply with Russian legislation the Instruction for use is supplemented by these Addendums, which are the parts of the instruction for use for the products supplied to the Russian Federation.

Addendum I to the document Instruction for use «KINEVO 900 Surgical Microscope with accessories»

Paragraph 2.2.1. "Intended use" stands for:

The KINEVO 900 is a surgical microscope intended for the illumination and magnification of the surgical area and for the support of visualization in surgical procedures.

KINEVO 900 is appropriate for use in cranial and spine neurosurgery, as well as otorhinolaryngology, in the field of auditory nerves and the base of the skull. In addition,

it is used for reconstructive and plastic procedures in traumatology, reconstructive and plastic surgery, oral surgery, chin surgery, and facial surgery in hospitals, clinics, and other medical institutions. KINEVO 900 is appropriate for interdisciplinary applications in microsurgery. When connected to the QEVO special camera, it is also intended for viewing the open operative field during general surgical procedures as well as for visual representation of the ventricles and structures inside the brain during surgical procedures in neurology. The use of INFRARED 800 with FLOW 800 (optional) during surgery applying fluorescent

technologies helps to view and visually assess bloodstream during surgery as well as vascular permeability in neurosurgical procedures, plastic and reconstructive surgery, as well as in coronary bypass surgery.

CAUTION! Injury to the patient's eye!

The device must not be used for ophthalmological procedures.

- Make sure that no xenon light and no laser radiation enters the patient's eyes.

Optional functions

The BLUE 400 option is designed for a surgical microscope and can be used in tumor surgery to visualize malignant tissue marked with fluorescence dyes. The YELLOW 560 option is an accessory for the intraoperative fluorescence excitation with light in the wavelength range from 460 nm to 500 nm. Light in the wavelength range between 540 and 690 nm can be observed.

Paragraph 2.2.2. "Indications " stands for:

2.2.2 Indications for use

Information on the indication of use and the target patient group is set out in the following table

Specialty:	• Neurosurgery and spine surgery • Orthopaedic surgery • Oral and maxillofacial microsurgery • P&R procedure • Not applicable for ophthalmology
Indications:	The ZEISS KINEVO 900 is a surgical

	microscope intended to illuminate and magnify the surgical area and to support visualization in surgical procedures.
Target patient population	No restrictions
Age	No restrictions
Sex	Both

Indications for use of BLUE 400, INFRARED 800 with FLOW 800 option and YELLOW 560 fluorescent modules

-The BLUE 400 option is designed for a surgical microscope and can be used in tumor surgery to visualize marked tumor tissue with fluorescence dyes.

-INFRARED 800 with FLOW 800 Option used during fluorescence guided surgery aids in viewing the visual assessment of intraoperative blood flow as well as vessel patency in surgical procedures in neurosurgery, plastic and reconstructive procedures and coronary artery bypass graft surgery.

-The YELLOW 560 option is an accessory for the intraoperative illumination of the surgical area with light in the wavelength range from 460 nm to 500 nm. Light in the wavelength range between 540 and 690 nm is particularly highlighted.

Contraindications and Possible Side Effects

Contraindications

KINEVO 900 Surgical Microscope

provided that ZEISS KINEVO 900 is used within its indication of use, there are no known contraindications.

Special endoscopic camera QEVO and Special endoscopic camera control unit QEVO ECU

The product is contraindicated for treatments in which the endoscopic method is contraindicated.

BLUE 400, YELLOW 560, INFRARED 800 with FLOW 800, fluorescent modules

provided that ZEISS modules are used within their indication of use, there are no known contraindications.

Side Effects

A high light intensity over a longer period of time can damage tissue and cause burns in the area of the illuminated surgical field.


Paragraph 3.3.2. " Approval label / UDI label" stands for:



Pos.	Explanation
1	Approval label

	■ Manufacturer: Carl Zeiss Meditec AG ■ Manufacturer's address: Goeschwitzer Strasse 51 - 52 07745 Jena, Germany ■ Manufacturer symbol:  ■ Device name: KINEVO 900 ■ Serial No. SN ■ Reference number REF ■ Date of manufacture, year-month-day ■ Date of manufacture symbol  ■ Rated voltage: 100 - 240V ■ Max. connected load: Max. 1350VA ■ Line frequency range: 50Hz - 60 Hz ■ Designation of origin for industrial products: MADE IN GERMANY ■ WEEE mark:  ■ CSA mark:  ■ CE marking ■ Maximum total mass: Max. 395 kg
2	UDI label ■ Machine-readable label (barcode) ■ UDI Device Identifier (UDI-DI) ■ UDI Production Identifier (UDI-PI)


Paragraph 3.3.5. " Labeling on stand, Part 3" stands for:


Pos.	Symbol	Explanation
1		Risk of crushing!
2		Approval label for INFRARED 800 (option) provides information on: ■ Manufacturer ■ Catalog number ■ Serial No. ■ Date of manufacture ■ License data ■ UDI marking: ■ CE marking
3		UDI label for ZEISS KINEVO 900, provides information on: ■ UDI marking ■ Serial No. ■ Device name
4		Fuse display F1
5		Read and observe the Instructions for Use.
6		Non-ionizing electromagnetic radiation
7		Potential equalization
8		Approval label for ZEISS KINEVO 900 provides information on: ■ Manufacturer ■ Catalog number ■ Serial No. ■ Date of manufacture ■ System data ■ Power data ■ UDI marking:

		■ Type of disposal ■ Total weight ■ CSA label ■ CE marking
9		SIP label, provides information on: ■ Manufacturer ■ Manufacturer's contact information ■ SIP number of device (SIP = Service Identification Program)
10	664010xxxx	Serial number of the device
11		MD labeling (Medical Device)* * Depending on the respective country- specific registration status
12		Port for navigation system
13		Identification label, options
		Identification for installed options

Paragraph 3.3.7. " Labeling on Light Source" stands for:

Pos.	Symbol	Explanation
1		Lamp replacement/manual lamp replacement (removal of light source cover)
2		Risk of crushing!
3		Hot surface
4		General warning
5		"Observe disposal regulations" label Do not dispose of electrical or electronic devices along with normal domestic waste.
6		Manual lamp change
7		Identification label, lamp 300 W+ xenon Contains the following information: ■ Manufacturer ■ Serial No. ■ Device name

		■ Catalog number
8		Light 1
9		Light 2
10		Identification label, "HLQ 300 Future" light source The identification label contains the following information: ■ Manufacturer ■ Serial No. ■ Device name ■ Catalog number

Paragraph 3.6.23 "Fluorescence" should be supplemented with the following information:

The fluorescent substances required for operation are not included in the scope of product delivery, but are supplied separately by other companies. Responsibility for quality, efficiency and safety of the fluorescent substances, state registration and compatibility with fluorescence modules BLUE 400, INFRARED with FLOW 800 option and YELLOW 560 lies on the manufacturer of the corresponding substances.

Only use the fluorescent substances that are registered in RF, for example, Aminolevulinic acid (registration certificate № ЛП-001848 dated 12.12.2018) for BLUE 400 or Fluorescein (registration certificate № ЛСР-000829/08 dated 28.01.2019) for YELLOW 560.

Paragraph 4.7. "Connecting the device to the power supply" stands for:

Connecting the protective ground!


CAUTION! To avoid the risk of electric shock, the device must only be connected to mains power that has a protective earth connection (IEC 60601-1)

► Connect the device only to electrical outlets that are properly equipped with a protective earth connection

Paragraph 5.1 "Safety during preparation" sub-paragraph "Do not cover the ventilation openings!" should be supplemented with the following information:

Ventilation openings

Paragraph 5.6 «Drape position» should be supplemented with the following information:

Drape position – factory (or arbitrarily changed by the user) setting, that allow to set the microscope in a comfortable position for attaching the ZEISS SMARTDRAPE sterile disposable drapes on the microscope (drape attaching), during daily startup.

Paragraph 7.2 «Contamination of the system» should be supplemented with the following information:

Special covers, supplied with the packaging, needed to close openings that are not in use.

Paragraph 10.1. "Regulatory information" sub-paragraph "Classification of the device in accordance with IEC 60601" stands for:

KINEVO 900 device Classification According to IEC 60601

- protection class against electric shock: 1 class
- operation mode: continuous.
- ingress protection (against water or solid particles): IP2X
- inappropriate for use in environments with high oxygen content

Paragraph 10.4. "Electrical data" replace the name

"Current consumption" with "operating power consumption"

"Automatic circuit breaker" with "Automatic fuse"

Paragraph 10.6 "Aiming beam laser" should be supplemented with the following Information:

- Beam divergence	beams are collimated	Not applicable
- Pulse pattern (pulse duration, repetition rate,)	radiation is continuous	Not applicable
- Maximum power or energy output	< 1 mW	

Laser beam function

Focus point is determined with two beam types, converging on the focal plane.

Paragraph 10.14.1 "Surgical microscope"

Instead of "Varioskop", read "working distance»

Paragraph 10.19.1.1 "Electromagnetic interference" supplemented with the following information:

See the test certificate in Appendix 3

Paragraph 12 "Disposal" supplemented with the following information:

Disposal in accordance with the sanitary rules and regulations of the Russian Federation. Medical waste hazard class: A.

For products ZEISS SMARTDRAPE sterile disposable drapes and VisionGuard replacement lenses:

Disposal in accordance with the sanitary rules and regulations of the Russian Federation. Medical waste hazard class: A.

Add "Software" to the instruction for use:

KINEVO 900 uses the current software version of 1.8 and higher.

Add "Warranty" to the instruction for use:

The warranty period is one year.

The warranty period of storage is 8 years.

The manufacturer disclaims all liability for claims arising from improper use of the product caused by non-compliance with the Instructions for Use.

Appendices:

Annex 2 – Composition list, accessories

Annex 3 - Main characteristics of components and accessories

Annex 4 - Labeling

Annex 5 -CB test Certificate

Annex 6 - Certificate of conformity

Composition list, accessories**KINEVO 900 Surgical Microscope with accessories:****1. KINEVO 900 Surgical Microscope, consisting of:**

- 1.1. Main unit of the KINEVO 900 surgical microscope (including optical microscope unit, microscope stand) with built-in unit for installation of fluorescent modules.
- 1.2. Wide-angle eyepiece (no more than 8 pcs. - if necessary), in variants:
 - 12.5x;
 - 10x.
- 1.3. Lamp unit (not more than 10 pcs. - if necessary).
- 1.4. Binocular tube (no more than 4 pcs. - if necessary), in variants:
 - 180° Tilttable tube;
 - Straight tube;
 - Foldable tube;
 - Foldable tube for mouth switch
- 1.5. ZEISS SMARTDRAPE sterile disposable drapes (in packing of 5 pcs. - not more than 80 packages, if necessary).
- 1.6. VisionGuard replacement lenses (in packing of 10 pcs. - not more than 20 packages, if necessary).
- 1.7. Special video camera, in variants:
 - Special video camera (302584-9721-000)
 - Special video camera (302584-9720-000)
 - Special video camera (302584-9772-000)
 - Special video camera (302584-9782-000)
- 1.8. Special video recorder (if necessary).
- 1.9. Power cord (no more than 5 pcs.).
- 1.10. Connection cables (no more than 25 pcs.), in variants:
 - Quad-SDI video cable, 10 m;
 - DVI-D cable, 5 m;
 - DVI-D cable, 10 m;
 - DVI-D cable, 0.3 m;
 - Displayport-D cable, 5 m;
 - Video Out cable, Lemo-E, 5 m;
 - Endo In cable, Lemo-E, 5 m;
 - Endo In cable, Lemo-DVI, 5 m;
 - Video Out cable (4K), ODU, 5 m;
 - HD-SDI cable, 10 m;
 - HDMI Video cable, 5 m.
- 1.11. Network adapter (if necessary).
- 1.12. Medium with special software.
- 1.13. Set of documentation, including:
 - Instruction for use «KINEVO 900 Surgical Microscope with accessories» with application.
 - Instruction for use «Special endoscopic camera QEVO and Special endoscopic camera control unit QEVO ECU » (if necessary)
 - QEVO Processing instruction (if necessary)
 - Instruction for use «Mouthswitch for 180° Tilttable tube» (if necessary).
 - Instruction for use «Multifunctional Foot Control Panel User Manual. Wireless FCP WL; wired FCP»(if necessary)

- Instruction for use «Fluorescent module INFRARED 800 with FLOW 800 option. For KINEVO 900» (if necessary).
- 1.14. Special endoscopic camera QEVO (no more than 5 pcs. - if necessary).
- 1.15. Special endoscopic camera control unit QEVO ECU (if necessary).
- 1.16. QEVO sterilization kit (no more than 5 pcs. - if necessary).
- 1.17. Second monitor (if necessary), in variants:
 - HD 2D monitor (305992-9061-000);
 - HD 3D monitor (305992-9060-000).
- 1.18. Rocker foot switch (if necessary).
- 1.19. Multifunctional foot control panel (if necessary), in variants:
 - wireless FCP WL;
 - wired FCP.
- 1.20. Connection cable for multifunctional foot control panel, 6m (if necessary).
- 1.21. Mouthswitch for 180° Tilttable tube (if necessary).
- 1.22. Stereo module for co-observation.
- 1.23. 3D glasses (no more than 50 pcs. - if necessary).
- 1.24. USB drive, 64 GB (no more than 6 pcs.).
- 1.25. External HDD, 1 TB.
- 1.26. Fluorescent module (no more than 4 pcs. if necessary), in variants:
 - 2.26.1 BLUE 400, including:
 - BLUE 400 target;
 - Filter in the observation path of the microscope for BLUE 400 (2 pcs.)
 - Excitation filter;
 - Fasteners;
 - Information sticker (2 pcs.);
 - 2.26.2 INFRARED 800 with FLOW 800 option, including:
 - Fluorescent INFRARED 800 target with FLOW 800 option;
 - Filter in the observation path of the microscope for INFRARED 800 with FLOW 800 (2pcs.);
 - Excitation filter;
 - Information sticker (2pcs.);
 - Fasteners;
 - 2.26.3 YELLOW 560, including:
 - YELLOW 560 target;
 - Filter in the observation path of the microscope for YELLOW 560 (2 pcs.)
 - Excitation filter;
 - Fasteners;
 - Information sticker (2 pcs.);
- 1.27. Fluorescent INFRARED 800 target with FLOW 800 option (no more than 2 pcs. - if necessary).
- 1.28. YELLOW 560 target (no more than 2 pcs. - if necessary).
- 1.29. BLUE 400 target (no more than 2 pcs. - if necessary).
- 1.30. License key for connecting navigation systems (if necessary).
- 1.31. License key for activating fluorescent module (no more than 5 pcs. - if necessary).
- 1.32. License key for activating DICOM connections (if necessary).
- 1.33. License key for wireless connection (if necessary).
- 1.34. License key for shared network


Accessories:

1. Adapter for external digital foto cameras (no more than 10 pcs.).
2. Dust cover (no more than 3 pcs.).

3. Optical cleaning kit (no more than 5 pcs.), consisting of: optical cleaning agent, dust removal brush, microfiber cloth, wet wipes (10 pcs.), kit bag.
4. Three-step magnification changer.

Main Specifications of Components and Accessories

Component/ Accessory	Weight (kg) ±5%	Dimensions ±5%	Figure	Other specifications																				
Main unit of the KINEVO 900 surgical microscope (including optical microscope unit, microscope stand) with built-in unit for installation of fluorescent modules	Max. 395 – Total weight of device with max. add. load	See the drawing Max. dimensions in transport position (W×H×D) 850×1900×1700 mm – microscope on a stand		Optical data: <table><tr><th>Description</th><th>Meaning</th></tr><tr><td>Working distance</td><td>About 200–625 mm</td></tr><tr><td>Zoom system magnification factor</td><td>0.4×–2.4×</td></tr><tr><td colspan="2">Laser with an aiming beam (autofocus feature)</td></tr><tr><td>Max. power (IEC 60825-1)</td><td><1 mW</td></tr><tr><td>Wavelength</td><td>λ 635–645 nm</td></tr></table> Additional Information ■ Magnification adjustment: motorized/manual ■ Focusing: motorized/manual ■ Sensor for working distance and magnification for neuronavigation ■ Motor-driven, movement speed adjustable, automatic adaptation to magnification Classification <table><tr><th>Description</th><th>Meaning</th></tr><tr><td>Liquid ingress and solid particle protection</td><td>IP2X</td></tr><tr><td>Operability in environments with high oxygen content</td><td>Not suitable for use in environments with high oxygen content</td></tr><tr><td>Operating mode</td><td>continuous operation</td></tr></table>	Description	Meaning	Working distance	About 200–625 mm	Zoom system magnification factor	0.4×–2.4×	Laser with an aiming beam (autofocus feature)		Max. power (IEC 60825-1)	<1 mW	Wavelength	λ 635–645 nm	Description	Meaning	Liquid ingress and solid particle protection	IP2X	Operability in environments with high oxygen content	Not suitable for use in environments with high oxygen content	Operating mode	continuous operation
Description	Meaning																							
Working distance	About 200–625 mm																							
Zoom system magnification factor	0.4×–2.4×																							
Laser with an aiming beam (autofocus feature)																								
Max. power (IEC 60825-1)	<1 mW																							
Wavelength	λ 635–645 nm																							
Description	Meaning																							
Liquid ingress and solid particle protection	IP2X																							
Operability in environments with high oxygen content	Not suitable for use in environments with high oxygen content																							
Operating mode	continuous operation																							

Electrical data

Description	Meaning
Rated voltage	100–240 V AC
Nominal power consumption	Max. 1350 V·A
Rated frequency	50–60 Hz
Protection class (as per EN 60601-1)	I
Remote connector	24 V
Fuse	Automatic fuse (Max. voltage 250V AC Rated current 16A) the higher the overcurrent, the shorter the break-time. (Also depends on the ambient temperature)

Laser with an aiming beam

Description	Meaning
Max. power (IEC 60825-1)	<1 mW
Wavelength	λ 635–645 nm
Laser class	II

Light source

Description	Meaning
Color temperature	5000 K (± 500 K)
Voltage	100-240V
Rated power	Approx. 300 W

Integrated HD camera

Description	Meaning
-------------	---------

Floor stand

Resolution	1920×1080 p
Signal-to-noise ratio	54 dB
Sampling frequency	50 Hz 59.94 frames/second

Integrated 4K camera

Description	Meaning
Resolution	3840×2160 p
Signal-to-noise ratio	54 dB
Sampling frequency	50 Hz 59.94 frames/second

Digital video outputs

Description	Meaning
Output: 6 HDMI/DVI	1920×1080 p 50/60
Output: 7, 8 Display port	1920×1080 p 50/60
Output: 9, 10 HD-SDI (-0.8 Vp-p / 75 ΩPAL)	1920×1080 p 50/60 1920x1080 i 50/60

Digital video inputs

Description	Meaning
Input: 3, 4 Display port	1920×1080 p 50/60
Input: 5 HDM/DVI	1920×1080 p 50/60 1280×720 p 50/60

4K video outputs

Description	Meaning
Outputs: 1, 2 HDMI UHD 4K	3840×2160 p 50/60

				<table><tr><td>Color sampling</td><td>4:2:0</td></tr></table> Mechanical data <table><tr><th>Description</th><th>Meaning</th></tr><tr><td>Max. additional load on microscope body</td><td>Max. 6 kg</td></tr><tr><td>Max. dimensions in transport position (W×H×D)</td><td>850×1900×1700 mm ±5%</td></tr><tr><td>Total weight of device with max. add. load</td><td>Max. 395 kg</td></tr><tr><td>Total weight of device incl. single-use transport box</td><td>Max. 525 kg</td></tr></table> Microscope movement in XY Microscope mode Tilt around X axis: –25° to +135° Rotation around Y axis: –45° to +45° (+/–0.5°) Microscope movement in XY Stand mode Movement depends on the orientation of the microscope If the microscope is oriented vertically downwards, Y range is the circle around the 1 axis. X range is determined by kinematics and tilts of axes 2 and 3. In horizontal position of the microscope, Y–X range is determined by kinematics and tilts of axes 2 and 3. Microscope movement speed in XY Microscope mode	Color sampling	4:2:0	Description	Meaning	Max. additional load on microscope body	Max. 6 kg	Max. dimensions in transport position (W×H×D)	850×1900×1700 mm ±5%	Total weight of device with max. add. load	Max. 395 kg	Total weight of device incl. single-use transport box	Max. 525 kg
Color sampling	4:2:0															
Description	Meaning															
Max. additional load on microscope body	Max. 6 kg															
Max. dimensions in transport position (W×H×D)	850×1900×1700 mm ±5%															
Total weight of device with max. add. load	Max. 395 kg															
Total weight of device incl. single-use transport box	Max. 525 kg															


				Min. speed: ≤ 1 °/s Max. speed: ≥ 11 °/s Microscope movement speed in XY Stand mode Min. speed: ≤ 0.025 m/s Max. speed: ≥ 0.06 m/s Focal length range Min. 200 mm (-10 mm) Max. 625 mm (+50 mm) Field of view size Field of view depends on the used combination of tube and eyepieces, as well as on current magnification and working distance. 170 mm tube and 10× eyepiece provide the following range: Max. magnification (2.4), min. working distance (200 mm): FOV 15 mm (+/-2.5%) Min. magnification (0.4), max. working distance (625 mm): FOV 225 mm (+/-2.5%) Illumination of the surgical field Illumination intensity range is 5–100% (linearity $\pm 15\%$) When the light source is set to 100%, the illumination in the center of the operating field is ≥ 120 klx when the system is set to white light mode at a working distance of 500 mm, with a diameter of 120 mm. Light-field size
--	--	--	--	--


				The diameter of the illuminated area depends on the working distance. Minimum working distance of 200 mm, the range is: Min. 15 mm (± 2 mm) Max. 100 mm (± 5 mm) For the maximum working distance of 625 mm the range is: Min. 26 mm (± 4 mm) Max. 270 mm (± 15 mm) Speed of motorized focusing from 5% to 100% Zoom speed 1% Speed of setting an operating mode <75 s Maximum force required to move <100 N Full magnification tolerance $\pm 7.5\%$ Divergence of right and left optical systems axes: vertical, up to 15' horizontal: - convergence, up to 45' - divergence, up to 10' Displacement of focus planes when changing magnification: ± 1.2 mm Focus difference between right and left optical systems: <2.1 mm Resolution in the center of the field: > 46 line pairs / mm at working distance 200 mm / maximum magnification
--	--	--	--	--



				> 27 line pairs / mm at working distance 400 mm / maximum magnification > 19 line pairs / mm at working distance 625 mm / maximum magnification Image tilt difference between right and left optical systems: $\leq 2^\circ$ Height difference of exit pupils between right and left optical systems: 1.5 mm at 0 diopters on the diopter scale Difference between the magnifications of the right and left optical systems: no more than 1,5% Minimum range of interpupillary distances: from 53 mm to 75 mm. Eyepiece diopter scale calibration allowance: ± 0.25 at 0 diopters on the diopter scale Minimal top of the eye eyepiece adjustment: from -8 to +5 available functions: - Performing a freeze frame (photo) - Starting and stopping video recording - Low-resolution video recording - Auto balance - Illumination control during focusing - Automatic brightness adjustment - Auto focusing - Laser focus
--	--	--	--	--


				- Replaceable UV/IR protective filters - Manual point lock (rotation) - Automatic point lock (rotation) - Position memory - Audio functionality: microphone, audio input/output ports, loudspeaker current type, power consumption / remote port AUX port for external signal, max. 24 V / 0.5 A Material ABS+pc HP-1001 thermoplastic (LOTTE ADVANCED MATERIALS), EN AC-Al Si12Cu1(Fe)F alloy (Multi-Wing) CuZn40Pb2 F43 alloy (Wetring) Brief skin contact. Nonsterile										
Eyepiece	Eyepiece 10×	0.216 kg	Diameter 46 mm, length 59 mm	 An eyepiece is an optical user interface between the surgeon's eye and the system. The system allows adjusting diopters to compensate for the surgeon's refractive error. Optical data of the widefield eyepiece (magnification 10×) <table><tr><th>Description</th><th>Meaning</th></tr><tr><td>Focal length</td><td>25 mm</td></tr><tr><td>Field of view</td><td>20.6 mm</td></tr><tr><td>Distance of the exit pupil from the last lens</td><td>24–25.5 mm</td></tr><tr><td>Diopter adjustment range</td><td>+5/–8</td></tr></table> Material: PP5CWZ PVC (Gummiwerk Kraiburg GmbH) HZK50 optical glass (CDGM GLASS) Brief skin contact. Nonsterile Order number 305542-0000-000	Description	Meaning	Focal length	25 mm	Field of view	20.6 mm	Distance of the exit pupil from the last lens	24–25.5 mm	Diopter adjustment range	+5/–8
Description	Meaning													
Focal length	25 mm													
Field of view	20.6 mm													
Distance of the exit pupil from the last lens	24–25.5 mm													
Diopter adjustment range	+5/–8													

	Eyepiece 12.5×	0.216 kg	Diameter 46 mm, length 59 mm		An eyepiece is an optical user interface between the surgeon's eye and the system. The system allows adjusting diopters to compensate for the surgeon's refractive error. Optical data of the widefield eyepiece (magnification 12.5×) <table><tr><th>Description</th><th>Meaning</th></tr><tr><td>Focal length</td><td>20 mm</td></tr><tr><td>Field of view</td><td>17.6 mm</td></tr><tr><td>Distance of the exit pupil from the last lens</td><td>22–23.5 mm</td></tr><tr><td>Diopter adjustment range</td><td>+5/–8</td></tr></table> Material: PP5CWZ PVC (Gummiwerk Kraiburg GmbH) HZK50 optical glass (CDGM GLASS) Brief skin contact. Nonsterile Order number 305543-9901-000	Description	Meaning	Focal length	20 mm	Field of view	17.6 mm	Distance of the exit pupil from the last lens	22–23.5 mm	Diopter adjustment range	+5/–8
Description	Meaning														
Focal length	20 mm														
Field of view	17.6 mm														
Distance of the exit pupil from the last lens	22–23.5 mm														
Diopter adjustment range	+5/–8														
Lamp unit		840 g	125×90 mm		The lamp container includes a xenon lamp as a main element of the illumination system. brand / order number / manufacturer HLQ300W / 304949-9002-000/ Carl Zeiss Meditec AG • Ray divergence 6° (10% half angle) through 25.4 mm exit window Light source - second spare xenon lamp - automatic lamp replacement - thermal insulation and UV filters - fluorescent filters <table><tr><th>Description</th><th>Meaning</th></tr><tr><td>Color temperature</td><td>Around 5600 K</td></tr><tr><td>Voltage</td><td>24 V</td></tr><tr><td>Rated output power</td><td>Approx. 300 W</td></tr></table>	Description	Meaning	Color temperature	Around 5600 K	Voltage	24 V	Rated output power	Approx. 300 W		
Description	Meaning														
Color temperature	Around 5600 K														
Voltage	24 V														
Rated output power	Approx. 300 W														

Binocular tube	180° tiltable tube	0.86	0.8×12.8×14 cm		<table><tr><th>Description</th><th>Meaning</th></tr><tr><td>focal length</td><td>170 mm</td></tr><tr><td>tilt angles</td><td>180° ±18°</td></tr><tr><td>interpupillary distance</td><td>55–75 mm</td></tr></table> Material ABS+pc HP-1001 thermoplastic (LOTTE ADVANCED MATERIALS), EN AC-Al Si12Cu1(Fe)F alloy (Multi-Wing) CuZn40Pb2 F43 alloy (Wetering) Brief skin contact. Nonsterile The tube is a second part of the optical interface between the surgeon and the system. It provides an interface for eyepieces and allows to adjust interpupillary distance and the tilt to ensure ergonomic working position for the surgeon. 303791-0000-000 / Carl Zeiss Meditec AG	Description	Meaning	focal length	170 mm	tilt angles	180° ±18°	interpupillary distance	55–75 mm
Description	Meaning												
focal length	170 mm												
tilt angles	180° ±18°												
interpupillary distance	55–75 mm												
Straight tube		0.528			<table><tr><th>Description</th><th>Meaning</th></tr><tr><td>focal length</td><td>170 mm</td></tr><tr><td>interpupillary distance</td><td>55–75 mm</td></tr></table> Material ABS+pc HP-1001 thermoplastic (LOTTE ADVANCED MATERIALS), EN AC-Al Si12Cu1(Fe)F alloy (Multi-Wing) CuZn40Pb2 F43 alloy (Wetering) Brief skin contact. Nonsterile 303765-0000-000 / Carl Zeiss Meditec AG	Description	Meaning	focal length	170 mm	interpupillary distance	55–75 mm		
Description	Meaning												
focal length	170 mm												
interpupillary distance	55–75 mm												


				The tube is a second part of the optical interface between the surgeon and the system. It provides an interface for eyepieces and allows to adjust interpupillary distance and the tilt to ensure ergonomic working position for the surgeon.										
Foldable tube	0.94	Width: 128.4 mm / Depth and height depend on the extension of the folding tube		<table><tr><th>Description</th><th>Meaning</th></tr><tr><td>focal length</td><td>switchable between 170 mm and 260 mm</td></tr><tr><td>tilt angles</td><td>170° ±2°</td></tr><tr><td>interpupillary distance</td><td>according to the scale 55–75 mm / min. 53 mm ±1 mm</td></tr><tr><td>magnification adjustment values (PROMAG)</td><td>approximately 50% magnification at a focal length of 260 mm</td></tr></table> Material ABS+pc HP-1001 thermoplastic (LOTTE ADVANCED MATERIALS), EN AC-Al Si12Cu1(Fe)F alloy (Multi-Wing) CuZn40Pb2 F43 alloy (Wetring) Brief skin contact. Nonsterile 303771-9021-000 / Carl Zeiss Meditec AG The tube is a second part of the optical interface between the surgeon and the system. It provides an interface for eyepieces and allows to adjust interpupillary distance and	Description	Meaning	focal length	switchable between 170 mm and 260 mm	tilt angles	170° ±2°	interpupillary distance	according to the scale 55–75 mm / min. 53 mm ±1 mm	magnification adjustment values (PROMAG)	approximately 50% magnification at a focal length of 260 mm
Description	Meaning													
focal length	switchable between 170 mm and 260 mm													
tilt angles	170° ±2°													
interpupillary distance	according to the scale 55–75 mm / min. 53 mm ±1 mm													
magnification adjustment values (PROMAG)	approximately 50% magnification at a focal length of 260 mm													

				the tilt to ensure ergonomic working position for the surgeon.											
Foldable tube for mouth switch	0.94	Width: 128.4 mm / Depth and height depend on the extension of the folding tube		<table><tr><th>Description</th><th>Meaning</th></tr><tr><td>focal length</td><td>switchable between 170 mm and 260 mm</td></tr><tr><td>tilt angles</td><td>170° ±2°</td></tr><tr><td>interpupillary distance</td><td>according to the scale 55–75 mm / min. 53 mm ±1 mm</td></tr><tr><td>magnification adjustment values (PROMAG)</td><td>approximately 50% magnification at a focal length of 260 mm</td></tr></table>	Description	Meaning	focal length	switchable between 170 mm and 260 mm	tilt angles	170° ±2°	interpupillary distance	according to the scale 55–75 mm / min. 53 mm ±1 mm	magnification adjustment values (PROMAG)	approximately 50% magnification at a focal length of 260 mm	Material ABS+pc HP-1001 thermoplastic (LOTTE ADVANCED MATERIALS), EN AC-Al Si12Cu1(Fe)F alloy (Multi-Wing) CuZn40Pb2 F43 alloy (Wetering) Brief skin contact. Nonsterile 303771-9110-000 / Carl Zeiss Meditec AG
Description	Meaning														
focal length	switchable between 170 mm and 260 mm														
tilt angles	170° ±2°														
interpupillary distance	according to the scale 55–75 mm / min. 53 mm ±1 mm														
magnification adjustment values (PROMAG)	approximately 50% magnification at a focal length of 260 mm														

The tube is a second part of the optical interface between the surgeon and the system. It provides an interface for eyepieces and allows to adjust interpupillary distance and the tilt to ensure ergonomic working position for the surgeon.

Fluorescence module	N/A Integrated	N/A Integrated	The unit is integrated BLUE 400  1) Target Blue 400 2) Filter in the observation path of the microscope for Blue 400 3) Excitation filter 4) Mount 5) Information sticker (2 pcs); INFRARED 800 with FLOW 800 option	BLUE 400 The BLUE 400 option is designed for a surgical microscope and can be used in tumor surgery to visualize malignant tissue marked with fluorescence dyes. <i>30BL400P / Carl Zeiss Meditec AG</i> Target Excitation filter: 400–410 nm Filters in the observation path of the microscope (2 pcs): 620–710 nm Mount <u>Illumination intensity:</u> minimal intensity for light source set at +100% shall be at least 13 mW/cm² in a focal plane in the center of the field of view. YELLOW 560 The YELLOW 560 option is an accessory for the intraoperative fluorescence excitation with light in the wavelength range from 460 nm to 500 nm. Light in the wavelength range between 540 and 690 nm can be observed. <i>30YE560P / Carl Zeiss Meditec AG</i> Target Excitation filter: 460–500 nm Filters in the observation path of the microscope (2 pcs): 540–690 nm Mount <u>Illumination intensity:</u> minimal intensity for light source set at +100% shall be at least 35 mW/cm² in a focal plane in the center of the field of view. INFRARED 800 <i>30IR800A / Carl Zeiss Meditec AG</i>
---------------------	-------------------	-------------------	--	---

- 1) Target INFRARED 800 with FLOW 800 option
- 2) Filter in the observation path of the microscope for INFRARED 800 with FLOW 800
- 3) Excitation filter
- 4) Information sticker (2pcs.)
- 5) Mount

YELLOW 560

Target

Filter in the observation path of the microscope:
815-920 nm

Excitation filter:
400-780 nm

Mount
information stickers

- magnification adjustment: computer-controlled automatic magnification and manual adjustment
- repeat frequency adjustment: 25 Hz / 29.97 Hz
- AutoGain
- AutoDetection
- AutoZoom
- Multivision mode


Filter for observation by the doctor:

Excitation filter:
D= 28 mm, weight 5 g
+/-5%


Fluorescent target:
93x90x19 mm, weight 0.26 kg
+/-5%


INFRARED 800 with FLOW 800 option:

- fluorescence brightness settings: automatic brightness mode
- fluorescence brightness gain
- AutoCorrection
- video playback
- current image recording

			 1) Target YELLOW 560 2) Filter in the observation path of the microscope for YELLOW 560 3) Excitation filter 4) Mount 5) Information stickers (2 pcs);	• time (scaling on the time axis. Imaging can be limited to the period from the start to the finish of fluorescence application) • AutoGain • AutoDetection • AutoZoom • gain adjustment: gain is controlled automatically or manually Range: 4 ... 100% • repeat frequency adjustment: 25 Hz / 29.97 Hz INFRARED 800 with FLOW 800 (option) can be used during surgery with fluorescence technology to help observe and visually evaluate the bloodstream, as well as blood vessel patency during neurosurgery, plastic, and reconstruction surgery, and heart bypassing.
Video camera, special	N/A Integrated	N/A Integrated		Video camera, special captures video streams from the left and right cameras to allow the following system functions: - Live video in GUI - photo capturing and video recording
Video camera, special	Control unit 2 kg Head unit 40 g	Control unit 30×23×4 cm Head unit 6×3.5×3.5 cm		Video cameras, special capture magnified image of the left and right microscope channels. Video camera, special (302584-9721-000) <i>Carl Zeiss Meditec AG</i> HD resolution (1920×1080) Video camera package consisting of: 1× Head Unit (CHU) 1× Camera Control Unit (CCU) 1× 2D VCU (Video control unit) 1× connecting cable Video camera, special (302584-9720-000) <i>Carl Zeiss Meditec AG</i>

			HD resolution (1920×1080) Video camera package consisting of: 2× Head Unit (CHU) 2× Camera Control Unit (CCU) 1× 3D VCU (Video control unit) 2× connecting cable Video camera, special (302584-9772-000) <i>Carl Zeiss Meditec AG</i> 4K UHD resolution (3840×2160) Video camera package consisting of: 1× 4K Head Unit (CHU) 1× 4K Camera Control Unit (CCU) 1× 2D VCU (Video control unit) 1× connecting cable Video camera, special (302584-9782-000) <i>Carl Zeiss Meditec AG</i> 4K UHD resolution (3840×2160) Video camera package consisting of: 2× 4K Head Unit (CHU) 2× 4K Camera Control Unit (CCU) 1× 3D VCU (Video control unit) 2× connecting cable • <i>Brightness adjustment, exposure time:</i> this parameter is available in automatic mode; the brightness is set automatically, while the brightness level can be set by the user from 0% to 100%. Manual mode: 1/50 s, 1/60... 1/10000 respectively • <i>color saturation:</i> the user can set it from -10 to +10. • <i>black level offset:</i> the user can set it from -10 to +10 • <i>red and blue adjustment:</i> the user can make fine red/blue adjustments from -100 to +100
--	--	--	--


			• <i>peak/average value level</i>: the setting is automatic • <i>white balance</i>: white balance procedure is automatic, apply manually <table><tr><th>Description</th><th>Meaning</th></tr><tr><td>Signal-to-noise ratio</td><td>54 dB</td></tr><tr><td>Sampling frequency</td><td>50 Hz 59.94 frames/second</td></tr></table>	Description	Meaning	Signal-to-noise ratio	54 dB	Sampling frequency	50 Hz 59.94 frames/second
Description	Meaning								
Signal-to-noise ratio	54 dB								
Sampling frequency	50 Hz 59.94 frames/second								
Second monitor	7.4	593.3×363.3 mm	Front view  Rear view HD 2D monitor (305992-9061-000) Carl Zeiss Meditec AG 24HD2DMO integrated monitor (HD 2D): - diagonal: 598.09 mm - display type: WUXGA LCD-TFT, LED backlight with a touchscreen - Response time (Tr + Tf): typical 30ms - viewing angle: typical H (horizontal): ±89° V (vertical): ±89°, min. H: ±85° V: ±85° at contrast ratio = 10:1 - display colors: >16.7 million. - resolution: 1920×1080 p - frequency: 50/60 Hz - panel brightness: typical 400 cd/m2; min. 360 cd/m2 at 0° viewing angle - contrast: typical 1000:1; min. 700:1 at 0° viewing angle - power supply parameters (voltage, power consumption, current type): input voltage 24 V DC ±15%, power consumption – max. 35 W at 100% brightness - rotation angle range: not applicable (rotation is impossible) HD 3D monitor (305992-9060-000) Carl Zeiss Meditec AG 24HD3DMO integrated monitor (HD 3D):						

				- diagonal: 598.09 mm - display type: WUXGA LCD-TFT - response time (Tr + Tf): typical 30ms - 2D viewing angle: typical H (horizontal): $\pm 89^\circ$ V (vertical): $\pm 89^\circ$, min. H: $\pm 85^\circ$ V: $\pm 85^\circ$ at contrast ratio = 10:1 - 3D viewing angle: Min. $\pm 5^\circ$, typ. $\pm 9^\circ$ @ 700mm viewing distance with <7% crosstalk - display colors: >16.7 million. - resolution: 1920×1080 p - frequency: 50/60 Hz - Panel brightness: typical 240 cd/m² - contrast: typical 1000:1; min. 700:1 at 0° viewing angle - power supply parameters (voltage, power consumption, current type): input voltage 24 V DC $\pm 15\%$, power consumption – 35 W. - rotation angle range: not applicable (rotation is impossible)
ZEISS SMARTDRAPE sterile disposable drapes	0.1	1320×3400 mm		Sterile covers for microscope protection during surgery Disposable, sterile Sterilization method – ethylene oxide Material: PE-HD Hostalen GC7260 high-density polyethylene (Bretthauer), A15-109-6011 dye blue (Albis) Plexiglas 8N polymethylmethacrylate (Degussa Röhm) ANTISTAT Blend low-density polyethylene (Microtek) 3M 9889 acrylic glue (3M) 11SMn30+C metal inserts + ISO 1456 -Fe/Ni20b nickel plating (Korrozin) SMART paper (Smart Tec) 306028-0000-000 / Microtec Dominicana, S.A. Air removal from the drape:

					Is done with an automatic suction pump, integrated into main unit of the microscope (supply voltage 24 V DC, air volume 16.2 m³ / h). More detailed information is in 5.7 Instructions for use «KINEVO 900 Surgical Microscope with accessories» Report on Validation of Sterilization, Report on Validation of Packaging of Products Subject to Final Sterilization – see Annex D Shelf life – up to 5 years
Connection cables	Quad-SDI video cable, 10 m	3.2	10 m		302584-8785-000 10 m-long cable for connecting an external monitor with KINEVO 900 via Quad – SDI interface ÖVE K41a EN50525 AS 3191 IEC 60227 polybutylene terephthalate (PBTP) Jacket: Polyvinyl chloride (PVC) by ekontor GmbH Heilbronn
	DVI-D cable, 5 m	0.9	5 m		305989-8698-000 For connection of external video devices (e.g. displays) via DVI-D for live video from the microscope. Jacket: Polyvinyl chloride (PVC) by ekontor GmbH Heilbronn
	DVI-D cable, 10 m	1.6	10 m		305989-8667-000 For connection of external video devices (e.g. displays) via DVI-D for live video from the microscope. ÖVE K41a EN50525 AS 3191 IEC 60227 polybutylene terephthalate (PBTP) Jacket: Polyvinyl chloride (PVC) by ekontor GmbH Heilbronn

DVI-D cable, 0,30m	0.1	30 cm		For connection of external video devices (e.g. displays) via DVI-D for live video from the microscope. ÖVE K41a EN50525 AS 3191 IEC 60227 polybutylene terephthalate (PBTP) Jacket: Polyvinyl chloride (PVC) by ekontor GmbH Heilbronn
Displayport-D Cable	0.9	5 m		302584-8746-000 For connection of external video devices (e.g. displays) via DisplayPort for live video from the microscope. ÖVE K41a EN50525 AS 3191 IEC 60227 polybutylene terephthalate (PBTP) Jacket: Polyvinyl chloride (PVC) by ekontor GmbH Heilbronn
The Video Out cable, Lemo-E, 5 m	0.85	5 m		302584-8743-000 For connection of external video devices (e.g. displays) via HDMI for live video from the microscope. ÖVE K41a EN50525 AS 3191 IEC 60227 polybutylene terephthalate (PBTP) Jacket: Polyvinyl chloride (PVC) by ekontor GmbH Heilbronn
- Endo In cable, Lemo-E, 5 m;	0.85	5 m		302584-8743-000 For connection of external endoscopes via HDMI for displaying endo video on KINEVO monitor (PiP or full screen). ÖVE K41a EN50525 AS 3191 IEC 60227 polybutylene terephthalate (PBTP) Jacket: Polyvinyl chloride (PVC) by ekontor GmbH Heilbronn
- Endo In cable, Lemo-DVI-, 5 m	0.9	5 m		302584-8745-000 For connection of external endoscopes via HDMI for displaying endo video on KINEVO monitor (PiP or full screen). ÖVE K41a EN50525 AS 3191 IEC 60227 polybutylene terephthalate (PBTP) Jacket: Polyvinyl chloride (PVC) by ekontor GmbH Heilbronn

Video Out cable (4K), ODU-, 5 m	0.9	5 m		302584-8758-000 ÖVE K41a EN50525 AS 3191 IEC 60227 polybutylene terephthalate (PBTP) Jacket: Polyvinyl chloride (PVC) by ekontor GmbH Heilbronn
HD-SDI cable, 10 m	0.75	10 m		10 m-long cable for connecting an external monitor to KINEVO 900 via SDI interface ÖVE K41a EN50525 AS 3191 IEC 60227 polybutylene terephthalate (PBTP) Jacket: Polyvinyl chloride (PVC) by ekontor GmbH Heilbronn
HDMI Video cable, 5 m.	0.9	5 m		302584-8758-000 ÖVE K41a EN50525 AS 3191 IEC 60227 polybutylene terephthalate (PBTP) Jacket: Polyvinyl chloride (PVC) by ekontor GmbH Heilbronn
Power cord	0.7	6 m		000000-0594-821 C19, 250VAC, 16A, 50/60Hz ÖVE K41a EN50525 AS 3191 IEC 60227 polybutylene terephthalate (PBTP) Jacket: Polyvinyl chloride (PVC) by ekontor GmbH Heilbronn

Network adapter	0.1 kg	42×22×93 mm		The WLAN network adapter is used to provide the functionality of a WLAN HotSpot or to connect the KINEVO 900 to a WLAN network. Transmitting and receiving frequencies: 2412–2472 MHz Transmission power the maximum output power may vary from country to country depending on local regulations.* Output power • 802.11b: +16 dBm min. • 802.11g: +14 dBm min. • 802.11a: +14 dBm min. • 802.11n HT20 (2.4 GHz): +13 dBm min. • 802.11n HT40 (2.4 GHz): +13 dBm min.
Medium with special software	0.01	70×20×10 mm		32 GB memory stick containing a microscope disk image and a program for disk cloning and disk recovery: Clonezilla
License key for activating the fluorescent module	0.01	70×20×10 mm		Memory stick with a license key
License key for connecting navigation systems	0.01	70×20×10 mm		Memory stick with a license key
License key for shared network	0.01	70×20×10 mm		Memory stick with a license key
Special endoscopic camera QEVO	250 g	Shaft: Length 120 mm Diameter 3.6 mm Cable length: At least 5 m		When the special endoscopic camera QEVO is connected it is intended for viewing the open operative field during general surgical procedures as well as for visual representation of the ventricles and structures inside the brain during surgical procedures in neurology. Upper segment, Shaft: EN 10088 1.4301 Stainless steel (manufacturers: Luminaries MATTHEYSA, RESONETICS SA (Medelec)) Aiming device, Cap: EN 10088 1.4301/1.4307 Stainless steel (manufacturers: Bebitz, ROLDAN SA) Flat glass: Sapphire glass (manufacturer: Badische Industrie-Edelstein-Gesellschaft)

1. The shaft contains rod lens optics. Light is conducted through rod lenses made of quartz glass and broken by air lenses between the rods.

1	Shaft	2	Hand grip
3	Cable	4	Connector

* The shaft contains rod lens optics. Light is conducted through rod lenses made of quartz glass and broken by air lenses between the rods.

The distal end of the shaft is equipped with an atraumatic tip and slanted 45° so that the surgical field can be viewed laterally.

** The hand grip contains a camera and a light source

Solder: AuSn 80/20 alloy (manufacturer: IPT-Albrecht GmbH)

Optical fiber: Quartz glass (Schott)

Glue: 2C-Epo-Tek OG214 epoxy resin (manufacturer: Epoxy Technology)

Handpiece: EN 10088 1.4307 Stainless steel (manufacturer: Mawatec Precision Mechanics)

Cable: Silicone (manufacturer: MediKabel)

Plug: Silicone. Holcosil LSR (manufacturer: ACC Silicones Limited, Holland Colours Europe B. V.)

Brief user skin and internal environment contacts.
Nonsterile, steam-sterilized

Parameters:

Rated voltage 24 V DC $\pm 10\%$

Connected load: max. 48 V·A

Rated connected load: 2 A

Luminous flux: 20–35 lm

Order number: 303155-9020-000, manufacturer: Carl Zeiss Meditec AG

Video Image Brightness Range GUI includes values from -5 to +5


operation cycle: 500 s / 1000 s

Direction of view at a distal end: 45°

resolution: 575 ± 50 TVL at a working distance of 10 mm


matrix type: sensor resolution = 1920×1080 pixels


field of view angle: 100° (-5°/+10°)

total length: 260 mm $\pm 5\%$


The product is sealed and resistant to biological fluids

				
Special endoscopic camera control unit QEVO ECU for Special endoscopic camera QEVO	2.2	240×272×59.3 mm		Control unit physically integrated in KINEVO 900 digital visualization node Only the front QEVO panel is available to the user Specifications: Rated voltage 24 V DC +/-10% Connected load: max. 48 V·A Rated connected load: 2 A Order number: 303155-9010-000, manufacturer: Carl Zeiss Meditec AG <u>Classification:</u> • Degree of protection against electric shock: Class I • Operability in environments with high oxygen content: inoperable in environments with high oxygen content • Operation mode: continuous operation <i>QEVO image angle alignment:</i> 8 different positions in a 360° circle in 45° increments

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1	Control unit	2	Socket for connecting the product plug
3	Connector panel	4	Control unit vent grill


white balance: The system features automatic white balance option

automatic switch-off function: the automatic switching function is available and can be activated in the GUI.

The light turns off after about 35 seconds if the QEVO camera is not in use

image brightness adjustment: -5...+5

marker function, image rotation settings: Current selected image orientation can be visualized in the video using the marker


Example:

**QEVO sterilization
box**

3.48

513×260×100

**304 stainless steel (Instru-Safe)
SP-65000 silicone (Summit medical)**


Brief user skin contact. Nonsterile, steam-sterilized

Multifunctional foot control panel

FCP WL multifunctional foot control panel, wireless

2.8 kg

410×243×128 mm

The multifunctional foot control panel is a user interface for controlling various KINEVO 900 functions with your foot.

304970-9200-000 / Carl Zeiss Meditec AG
(FCP WL multifunctional foot control panel, wireless)

Transmission frequencies from 2402 to 2480 MHz
Receiving power from -82 dBm to 0 dBm
Transmission power 0 mW – max. 2.5 mW – Class 2
Current consumption from 100 mA

HFC control panel (FCP WL)

power supply parameters (battery type, power supply voltage, number of batteries, operating time):

3 alkaline manganese acid type C batteries

Foot control panel is designed for continuous operation

IP classification: IPX8

304970-9100-000 / Carl Zeiss Meditec AG
(FCP multifunctional foot control panel, wired):

Has a 6 m-long cable

15 V DC $\pm 10\%$ via suspension system

Current consumption: max. 200 mA

power supply parameters (battery type, power supply voltage, number of batteries, operating time):


3 alkaline manganese acid type C batteries

Designed for continuous operation

IP classification: IPX8

14 wireless and wired functions:

FCP
multifunct
ional foot
control
panel,
wired

Longitudinal configuration functions	Transverse configuration functions
The joystick provides control in four following directions: 1 – forward 2 – backward 3 – left 4 – right 5 – Switch 6 – Switch 7 – Focus Down 8 – Focus Up 9 – Switch 10 – Switch 11 – Zoom In 12 – Zoom Out 13 – Switch 14 – Switch	The joystick provides control in four following directions: 1 – forward 2 – backward 3 – left 4 – right 5 – Switch 6 – Switch 7 – Focus Down 8 – Focus Up 9 – Switch 10 – Zoom Out 11 – Zoom In 12 – Switch 13 – Switch 14 – Switch

Multifunction foot
control panel
connection cable

0.2

6 m


304970-8760-000

Multifunction foot control panel connection cable

Mouthswitch for
180° Tilttable tube

0.47 kg

200×70×50 mm

Mouthswitch for 180° Tilttable tube is an additional simple user interface that allows to control the microscope with a mouth, hands-free. Pressing together the two parts of the mouth switch will release the magnetic brakes of the suspension system and enable you to guide the balanced surgical microscope almost effortlessly.

Material

1.4301 grade stainless steel as per EN 10088 (Haager GmbH & Co KG)

Thermoplastic ABS+PC HP-1001 (LOTTE ADVANCED MATERIALS)


Brief user skin contact. Nonsterile

Order number 000000-1177-805 / Carl Zeiss Meditec AG


connection cable length: 540 mm ±5%

height adjustment: 175–195 mm (depends on the tilt angle)

tilt angle range:

mouthpiece distance adjustment: limited by the interface surface at the tube, around 16 mm behind it at 85 mm in front of it for tilt mouth switch at 0°

approximately up to 115mm in front of the surface for tilt switch at 45°

Power supply parameters (voltage, current type, connected load/current): no power/current consumption
compression load of two parts: <2 H
Disable Magnetic Brakes of the Suspension System function

Care of the unit:


Cleaning mechanical surfaces



All mechanical surfaces of the equipment can be cleaned by wiping with a moist cloth. Do not use any aggressive or abrasive cleaning agents.

Wipe off any residue with a mixture of 50% ethyl alcohol and 50% distilled water plus a dash of household dish-washing liquid.


Place of installation is in 2.2.4.5. of this document (see position 9), installation of the mouthswitch is in 4.5 Instructions for use «KINEVO 900 Surgical Microscope with accessories»

Order of application is in Instructions for use «Mouthswitch for 180° Tilttable tube»


INFRARED 800 target	0.26 kg	93×90×19 mm		Target for function testing of the fluorescent module INFRARED 800 with FLOW 800 option 302581-9247-000 Order of application is in 5.1 Instructions for use «Fluorescent module INFRARED 800 with FLOW 800 option. For KINEVO 900»
License key for establishing a wireless connection	0.01 kg	70×20×10 mm		Memory stick with a license key
Stereo glasses	0.2	145×47×158		000000-1992-943 (5 glasses) PH38CP polycarbonate, manufacturer: Fa. Willbag Brief user skin contact. Nonsterile For viewing images on a 3D monitor
Rocker foot switch	0.752	190×79×29.5 mm, with a 3 m cable		Rocker foot switch is a user interface for controlling various KINEVO 900 functions with your foot. 2 functions (rocker foot switch has two buttons that can be activated (pressed) or deactivated (not pressed). KINEVO 900 software allows you to set these two buttons for various system functions. Order number 305989-8609-000 / Steute Medizintechnik Max. voltage 25 V AC / 60 V DC Max. current: 1 A Max. connected load 30 W IP classification: IPX8


				
Adapter for external digital foto cameras	0.42 kg	200×70×53 mm	CSA mark (Canadian Standards Association):  1/6164 	Adapter is intended for connecting external cameras to the surgical microscope. You can optionally connect a Canon DSLR camera to the device via a photo adapter and take photos with it. Photos are snapshots and can be taken in various situations: • During an operation when no video recording is running • During an operation when a video recording is running Note: The photos taken are saved exclusively to the storage medium of the DSLR camera. Place of installation of the adapter is in 2.2.4.5 of this document (see position 10). Material: AlMg5 Order of application is in 6.10 Instructions for use «KINEVO 900 Surgical Microscope with accessories»


				
BLUE 400 Target	0.26	93×90×19 mm		Fluorescent target for function testing of the fluorescent module BLUE 400 302581-9052-000 Order of application is in 6.14.2 Instructions for use «KINEVO 900 Surgical Microscope with accessories»

YELLOW 560 target	0.26	93×90×19 mm		Fluorescent target for function testing of the fluorescent module YELLOW 560 302582-9208-000 Order of application is in 6.15.2 Instructions for use «KINEVO 900 Surgical Microscope with accessories»
VisionGuard replacement lenses	0.555	5×5×2 cm		Plexiglas 8N polymethylmethacrylate (Degussa Röhm) Sterile Ethylene oxide sterilization. Report on Validation of Sterilization, Report on Validation of Packaging of Products Subject to Final Sterilization – see Annex D 306001-0000-000 / Carl Zeiss Meditec AG Shelf life – up to 5 years
License key for activating the DICOM connections	0.01	70×20×10 mm		Memory stick with a license key

Dust protection cover		215 g	800×1200 mm		To protect the microscope from dust
USB stick – 64 GB		0.01	70×20×10 mm		Size: 64 GB 000000-0576-392
1 TB external hard drive		0.26	110×80×15 mm		Size: 1 TB 000000-0574-996
Lens cleaning kit	Lens cleaning tool	0.16 kg	100×100×100 mm		Volume: 30 ml; Ingredients: diethylene glycol
(order number	Dust brush		l×d: 18×18 mm		-


000000-2096-685)	Microfiber cloth		18×18 cm		Material: microfiber
	Wet wipes (10 pcs.)		50×55 mm		Material: microfiber; impregnation agent composition: isopropanol, water
	Kit bag		100×110×40 mm		-
Three-stage zoom switch	0.448 kg	12×6×6 cm	 	303429-9903-000 Three-stage zoom switch is an optical accessory that can be installed between the microscope and the tube. It contains optics for additional magnification. Switch stages 0.6× 1× 1.6× Material: AlMg3	

				
Stereo module for co-observation	1084 g	Ø70×300 mm		000000-1063-869 Material ABS+pc HP-1001 thermoplastic (LOTTE ADVANCED MATERIALS), EN AC-Al Si12Cu1(Fe)F alloy (Multi-Wing) CuZn40Pb2 F43 alloy (Wetering) Brief skin contact. Nonsterile Purpose: co-observation (1) Connection (insert it into the beam splitter of the surgical microscope) (2) Axis for moving the stereo module (3) Release lever (4) Axis for moving the stereo module (5) Dovetail ring mount (6) Ring (7) Fastening screw (8) Axis for moving the stereo module

Packaging and product labeling

The KINEVO 900 system has the following labels on its shipping packaging and the devices.

Labeling on the Microscope

Pos.	Symbol	Explanation
1		Max. load: 6 kg
2		Open, right/left co-observation port
3		Warning of laser beam (labeling for autofocus option)
4	Varioskop	Brand name
5		Brand / ZEISS logo
6	Z	Zoom
7	F	Focus
8	302584-9000-000 XXXXXXXXXX	Microscope identification label (inside)
9		Switch over, right/left co-observation port
10		Luminous-field diaphragm
11		Focus/working distance
12		Adjust zoom magnification manually

Approval label / UDI (unique device identifier) label

Pos.	Explanation
1	Approval label ■ Manufacturer: Carl Zeiss Meditec AG ■ Manufacturer's address: Goeschwtitzer Strasse 51 - 52 07745 Jena, Germany ■ Manufacturer symbol: ■ Device name: KINEVO 900 ■ Serial No. SN ■ Reference number REF ■ Date of manufacture, year-month-day ■ Date of manufacture symbol ■ Rated voltage: 100 - 240V ■ Max. connected load: Max. 1350VA ■ Line frequency range: 50Hz - 60 Hz ■ Designation of origin for industrial products: MADE IN GERMANY ■ WEEE mark: ■ CSA mark: ■ CE marking ■ Maximum total mass: Max. 395 kg
2	UDI label ■ Machine-readable label (barcode) ■ UDI Device Identifier (UDI-DI) ■ UDI Production Identifier (UDI-PI)

Project of labeling on the territory of the Russian Federation

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
Labeling on stand, Part 1

Pos.	Symbol	Explanation
1	KINEVO 900	Device name
2		Brand / Logo
3		Seal for drape suction
4		Risk of crushing!
5		Transport position: • Observe the Instructions for Use • Pushing position • General warning


Labeling on stand, Part 2



Pos.	Symbol	Explanation
1		Standby/ON-OFF switch
2		RFID scanner or reading device
3	USB OPEN	Open USB cover
4		4K video outputs • Number 1: 4K HDMI output • Number 2: 4K HDMI output output (present only on 4K option)
5		Video input • Number 3: Display port input • Number 4: Display port input • Number 5: HDMI/DVI input
6		Video output • Number 6: HDMI/DVI output • Number 7: Display port output • Number 8: Display port output • Number 9: HD-SDI / 3G-SDI output • Number 10: HD-SDI / 3G-SDI output
7		Remote port
8		Foot control panel, 14 functions
9		Audio in
10		Audio out

11	LAN	Ethernet port (LAN)
12		Rocker switch
13		4K video output (present only on 4K option, not yet available) • Number 11: 3G-SDI output, left (2D/3D) • Number 12: 3G-SDI output, left (2D/3D) • Number 13: 3G-SDI output, left (2D) • Number 14: 3G-SDI output, left (2D) • Number 15: 3G-SDI output, right (3D) • Number 16: 3G-SDI output, right (3D)

Labeling on stand, Part 3

Pos.	Symbol	Explanation
1		Risk of crushing!
2		Approval label for INFRARED 800 (option) provides information on: ■ Manufacturer ■ Catalog number ■ Serial No. ■ Date of manufacture ■ License data ■ UDI marking: ■ CE marking
3		UDI label for ZEISS KINEVO 900, provides information on: ■ UDI marking ■ Serial No. ■ Device name
4		Fuse display F1

5		Read and observe the Instructions for Use.
6		Non-ionizing electromagnetic radiation
7		Potential equalization
8		Approval label for ZEISS KINEVO 900 provides information on: ■ Manufacturer ■ Catalog number ■ Serial No. ■ Date of manufacture ■ System data ■ Power data ■ UDI marking: ■ Type of disposal ■ Total weight ■ CSA label ■ CE marking
9		SIP label, provides information on: ■ Manufacturer ■ Manufacturer's contact information ■ SIP number of device (SIP = Service Identification Program)
10	664010xxxx	Serial number of the device
11		MD labeling (Medical Device)* * Depending on the respective country- specific registration status
12		Port for navigation system
13		Identification label, options
		Identification for installed options


Labeling on stand base




Pos.	Symbol	Explanation
------	--------	-------------

1		Pedal locking tabs
2		Pedal for straight-ahead movement

Labeling on Light Source


Pos.	Symbol	Explanation
1		Lamp replacement/manual lamp replacement (removal of light source cover)
2		Risk of crushing!
3		Hot surface
4		General warning
5		“Observe disposal regulations” label Do not dispose of electrical or electronic devices along with normal domestic waste.
6		Manual lamp change
7		Identification label, lamp 300 W+ xenon Contains the following information: ■ Manufacturer ■ Serial No. ■ Device name ■ Catalog number

8		Light 1
9		Light 2
10		Identification label, "HLQ 300 Future" light source The identification label contains the following information: ■ Manufacturer ■ Serial No. ■ Device name ■ Catalog number


Labeling for autofocus option

Pos.	Symbol	Explanation
1		Warning of laser beam
Country-specific labels: Laser class II – Laser radiation – Do not look directly into the beam or directly observe it with optical devices		
2	 	Country-specific labels: Labeling in accordance with IEC 60825-1:2007 and IEC 60825-1:2014 Laser radiation - Do not look into the beam or view directly with optical instruments Performance features • Laser class: II • Power output max: < 1 mW • Wavelength: λ 635 - 645
3		Label in Russian: Class 2 laser USA label: Laser class 2

Labeling on packaging


Symbol	Symbol	Explanation
--------	--------	-------------

	Indication of direction "This side up"	Indicates the correct upright position of the package.
	Fragile	Handle with care
	Keep dry	Protect packaging and packaged contents from wetness.
	Do not stack	Stacking of the packages is not permitted. No load should be placed on the package.
	Permissible temperature	The product may only be transported and stored at a temperature range of min. -20°C to max. +60°C.
	Packing unit	Specification of the number of packing units
	Permissible relative humidity	The product may only be transported and stored at humidity of min. 10% and max. 90% RH
	Permissible atmospheric pressure	The product may only be transported and stored at an atmospheric pressure of min. 500 hPa and max. 1060 hPa.
KINEVO 900		Device name
		Brand / Logo
		Min. 1.2 m forklift length
	Center of gravity	Indicates the package's center of gravity.
		Shockwatch
		Tiltwatch

Label on ZEISS SMARTDRAPE sterile disposable drapes

For transportation


Symbols on the label

	Do not re-use
	Sterilized using ethylene oxide
	Do not re-sterilize
	Does not contain latex
	Number of products
	Use-by date
	Date of manufacture
	Consult instructions for use
	Keep dry
	Keep away from sunlight


	Do not use if package is damaged
	Batch code
	CE label with the notified body's number
	Manufacturer


VisionGuard replacement lenses marking

Symbol	Explanation
	Catalog number
	Does not contain latex
	Keep away from sunlight
	Do not use if package is damaged
	Keep dry
	Batch code
	Sterilized using ethylene oxide
	Do not re-use
	Do not re-sterilize
	Use-by date
	Date of manufacture
	Consult instructions for use
VisionGuard Replacement Lenses	Device name
	Manufacturer
	Logo
	CE label with the notified body's number

QEVO camera labeling


1	INT 500s/1000s	Intermittent flow	Operating mode: do not use [▶ 27] the device in continuous mode.
2		UDI	UDI marking:
3		Catalog number	Reference number
4		CE marking	The medical device description corresponds to Annex XII of the Medical Devices Directive 93/42/EWG.
5		Type BF Applied Part	Meets BF requirements according to IEC 60601-1 for protection against electric shock.
6		Serial No.	
7		Maintenance intervals	Information about the frequency of maintenance.

Power plug labeling




Pos.	Symbol	Name	Explanation
1		Arrow	Power plug connection direction on the control panel of the QEVO ECU Terminal Box
2		Type BF Applied Part	Meets BF requirements according to IEC 60601-1 for protection against electric shock.


Labeling on packaging

1		Keep away from sunlight	Protect the packaging and its contents from high temperatures and sunlight.
2		Keep dry	Protect packaging and packaged contents from wetness.
3		Operating instructions	Read and follow the operating instructions.
4		Nonsterile	The device in the package is clean but not sterile. The product must be completely treated before first use.
5		-	The product is only allowed to be used by doctors or other persons upon a doctor's instructions.

Pos.	Symbol	Name	Explanation
6		MR unsafe	The device can not be stored or used in an environment where magnetic resonance occurs.
7		UDI	UDI marking:
8		Unique device identifier number	Contains information for unambiguous identification of the device.
9		CE marking	The medical device description corresponds to Annex XII of the Medical Devices Directive 93/42/EEG.
10		Endoscopic	Device type
11	Made in Germany	Manufacturer's country	The device is made in Germany.
12		Packing unit	The packaging contains one device
13		Atmospheric pressure limitation	The product may only be transported and stored at an atmospheric pressure of min. 500 hPa and max. 1060 hPa.
14		Air humidity limitation	The product may only be transported and stored at humidity of min. 10% and max. 92% RH.
15		Temperature limitation	The product may only be transported and stored at a temperature range of min. -20°C to max. +60°C.
16		Serial No.	-

Pos.	Symbol	Name	Explanation
17	QFVO	Device name	-
18		Catalog number	Reference number
19		Manufacturer	Contains the name and address of the manufacturer and the manufacture date of the device.
20	INT 500s/1000s	Intermittent duty	Operating mode: do not use [▶ 27] the device in continuous mode.


Pos.	Symbol	Name	Explanation
1		Top direction for transportation and storage of the packaging	The arrows should point up; the packaging is placed correctly here.
2		Fragile: transport carefully.	The contents of the packaging may break.

Pos.	Symbol	Name	Explanation
------	--------	------	-------------

3		Keep dry	Protect packaging and packaged contents from wetness.
4		Air humidity limitation	The product may only be transported and stored at humidity of min. 10% and max. 92% RH.
5		Atmospheric pressure limitation	The product may only be transported and stored at an atmospheric pressure of min. 500 hPa and max. 1060 hPa.
6		Temperature limitation	The product may only be transported and stored at a temperature range of min. -20°C to max. +60°C.
7		Number of packaging units	The packaging contains one device

Type plate with control unit data

1	Carl Zeiss Meditec AG	Manufacturer
---	-----------------------	--------------

2	Goeschitzer Strasse 51-52 07745 Jena, Germany	Manufacturer's address:
3	max. 48 VA	Power consumption
4	CE marking	The medical device description corresponds to Annex XII of the Medical Devices Directive 93/42/EEG.
5	CSA label	The medical device complies with the US and Canadian safety regulations for medical electrical systems.
6	Comply with disposal regulations	Do not dispose of electrical or electronic devices along with normal domestic waste. See the Disposal section for more detailed instructions on how to dispose of used electrical and electronic devices.
7	MADE IN GERMANY	Manufacturer's country
8	24 V DC	Rated voltage
9	QEVO ECU	Control unit name

Product identification plate

1	MD	Label for a medical device
2	UDI	UDI marking:
3	Unique device identifier number	Contains information about the manufacture date and the current serial number
4	Date of manufacture	Year-month-day
5	SN	Serial No.
6	REF	Catalog number

Multifunctional foot control panel labeling

1A

1 Type plate

Type plate contains the following information:

- Manufacturer symbol 
- Catalog number 
- Product name (wired version) FCP
- Manufacturer (company name) Carl Zeiss Meditec AG
- Manufacturer address Goeschwitzer Strasse 51–52
07745 Jena, Germany
- Serial No. 
- Product name (wireless version) FCP WL
- Protection rating IPX8

1B


2 Follow the Requirements of the Manual label

Follow the manual or other applicable documents.

3 Follow the Rules of Disposal label

Do not dispose of electrical and electronic devices as normal municipal waste. For more information on disposal of electrical and electronic devices, see the Maintenance section.

4 Date of manufacture

This symbol indicates the product's date of manufacture.

5 Label confirming the radio frequency approval

See the document G-30-2021 (Information on Radio Frequency Approval).

6 Radio frequency approval label at the base of the suspension system

See the document G-30-2021 (Information on Radio Frequency Approval).

7 Do Not Use Rechargeable Batteries label

8 Required Number and Type of Batteries label

3 × LR 14 (C battery)


9 Battery Orientation, Correct Polarity label


10 Pairing label

A label on the suspension system indicating the wireless connection of the foot control panel to the suspension system.


180° Tilttable Tube label


Straight Tube label


Label

Foldable tube


Foldable tube for mouth switch

Symbols	Explanation
 	(1) Direction of rotation
	(2) Zoom adjustment knob
Carl Zeiss f 170 / f 260	(3) Information about the manufacturer and the product
Carl Zeiss Meditec AG f 170 / f 260 303771-9XXX-000	(4) Label containing: - Manufacturer - Product information (tube focal length) - Order number
303771-9XXX SN XXXXXX 	(5) Barcode label

Fluorescent target INFRARED 800 with FLOW 800 option label

BLUE 400 Target label


YELLOW 560 Target label

Explanation


Symbols	Explanation
Carl Zeiss Meditec AG	Information about the manufacturer
Target INFRARED 800	Information about the product
302581-9247-000	Catalog number

Monitor Mobile Stand label


Symbol	Explanation
Brand name: pro-cart	Information about the product
Article no.: KU.5002.800	Catalog number
Serial / Lot no.: 001-091812	Series/lot number
ITD	Manufacturer's logo
CE	CE marking

Three-stage Zoom Switch label

Symbols	Explanation
0.6/1/1.6	Zoom indication

Stereo Module for Co-Observation label

Symbols	Explanation
	Information about the manufacturer Catalog number CE marking


Rocker Foot Switch label

Symbol	Explanation
--------	-------------

	Serial No.
	CSA symbol (for Canada and the United States)
	CE marking

Mouthswitch for 180° Tiltable tube label

Label contents:
- Manufacturer
- Order number
CE marking

Ref. Certif. No.

CA/21296/CSA

IEC SYSTEM FOR MUTUAL RECOGNITION OF TEST CERTIFICATES FOR ELECTRICAL EQUIPMENT (IECEE) CB SCHEME

CB TEST CERTIFICATE

Product

Surgical microscope on floor stand

Name and address of the applicant

Carl Zeiss Meditec AG
Goeschwtitzer Strasse 51-52, 07745 Jena, Germany

Name and address of the manufacturer

Carl Zeiss Meditec AG
Goeschwtitzer Strasse 51-52, 07745 Jena, Germany

Name and address of the factory

Carl Zeiss Meditec AG
Rudolf-Eber-Strasse 11, 73447 Oberkochen, Germany

Note: When more than one factory, please report on page 2

☐ Additional Information on page 2

Ratings and principal characteristics

KINEVO 900: 100 – 240 V, 50 – 60 Hz, max. 1350 VA
(incl. INFRARED 800)

QEVO ECU: 24 V DC, max. 48 VA (incl. QEVO)

Trademark (if any)

Customer's Testing Facility (CTF) Stage used

CTF Stage 3

Model / Type Ref.

KINEVO 900 (REF-No. 6640)
INFRARED 800 (REF-No. 7012)
QEVO ECU (REF-No. 7213)
QEVO (REF-No. 7212)

Additional information (if necessary may also be reported on page 2)

Additional test of Immunity to electromagnetic disturbances created by HF surgical equipment acc. to IEC 60601-2-2:2017, Annex BB.4. INFRARED 800, QEVO and QUEVO ECO are parts / built in options of the KINEVO 900 MEE.

☐ Additional Information on page 2

A sample of the product was tested and found to be in conformity with

IEC 60601-1-2: 2014

As shown in the Test Report Ref. No. which forms part of this Certificate

CB 70176409 + Enclosure 1

This CB Test Certificate is issued by the National Certification Body

CSA
Group

CSA Group
178 Rexdale Boulevard
Toronto, ON M9W 1R3 Canada

Signature: Michel Brossoit, P.Eng.

Date: 2018-04-11

Certificate of Compliance

Certificate: 70121971

Master Contract: 176164

Project: 80013302

Date Issued: 2019-10-08

Issued To: Carl Zeiss AG
22 Carl-Zeiss-Strasse
Oberkochen, Baden-Württemberg, 73447
Germany

Attention: Axel Keim

The products listed below are eligible to bear the CSA Mark shown with adjacent indicators 'C' and 'US' for Canada and US or with adjacent indicator 'US' for US only or without either indicator for Canada only.

Issued by: Ronald Garcia
Ronald Garcia

PRODUCTS

CLASS - C878001 - MEDICAL ELECTRICAL EQUIPMENT/SYSTEMS

CLASS - C878081 - MEDICAL ELECTRICAL EQUIPMENT/SYSTEMS Certified to US Standards

Medical Electrical Equipment, Model/Type: KINEVO 900, REF 6640, Appliance coupler, Mobile, rated:
100-240 Vac, 50-60 Hz, max. 1350 VA

1. Medical device protection against electric shock: Class I
2. Applied Part protection against electric shock: No applied part; with option "QEVO" Type BF
3. Degree of protection against ingress of water or particulate matter: No degree of protection
4. Method of Sterilization: None
5. Suitability for use in an Oxygen Rich Environment: Medical device not intended to be used in an Oxygen Rich Environment

Certificate: 70121971

Project: 80013302

Master Contract: 176164

Date Issued: 2019-10-08

6. Suitability to use Medical device in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide: Medical device not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide.
7. Mode of operation: Continuous
8. Environmental Conditions: Normal: 10-40°C, 30-75% RH, 700-1060hPa
9. Software Version : 1.6

APPLICABLE REQUIREMENTS

CSA Standards:

CAN/CSA-C22.2 No. 60601-1:14	CAN/CSA-C22.2 No. 60601-1:14: Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance (Adopted IEC 60601-1:2005 edition 3.0 + AMENDEMENT 1, 2012-07, MOD)
CAN/CSA-C22.2 NO. 60601-1-6:11 + AMD1	Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral standard: Usability (Adopted IEC 60601-1-6:2010, third edition, 2010-01 + A1:2013 – edition 3.1 This consolidated version consists of the third edition (2010) and its amendment 1 (2013-10)).
CAN/CSA-E60825-1:15	Safety of laser products — Part 1: Equipment classification and requirements (Adopted IEC 60825-1:2014, third edition, 2014-05)
CAN/CSA-C22.2 No. 62471:12	Photobiological safety of lamps and lamp systems (Adopted IEC 62471:2006, first edition, 2006-07)

ANSI/AAMI Standards:

ANSI/AAMI ES60601-1:2005/(R)2012, AND C1:2009 AND A2:2010(R)2012 (Consolidated text - edition 3.1)	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance (IEC 60601-1:2005+A1:2012, MOD).
IEC 60601-1-6:2006	Medical electrical equipment, Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability
IEC 60825-1:2014	Safety of laser products - Part 1: Equipment classification and requirements
IEC 62471:2006	Photobiological safety of lamps and lamp systems

Subject to the following qualifications:

- (1) This report describes the certification of the Medical Electrical Equipment or the Medical Electrical System with a North American Certified power supply cord set as indicated in the CSA description report.
- (2) A Medical Electrical Equipment or Medical Electrical System for shipping to outside of North America and not provided with a North American power supply cord set as described in the report, is certified as a component or a sub-assembly, and this certification does not cover the power supply cord set.
- (3) The user replaceable mains (line) fuse must be an approved type acceptable to the authorities where the equipment is sold.

Certificate: 70121971
Project: 80013302

Master Contract: 176164
Date Issued: 2019-10-08

- (4) Evaluated to CAN/CSA-C22.2 No. 60601-1:14 and ANSI/AAMI ES60601-1:2005/(R)2012, AND c1:2009 AND a2:2010(r)2012 (Consolidated text - edition 3.1) excluding (Not Evaluated) requirements for Electromagnetic compatibility (Clause 17), Biocompatibility (Clause 11.7)
- (5) SAFETY HAZARDS resulting from the intended physiological function of EQUIPMENT covered by this Standard are not considered.
- (6) Interconnection of this medical device with other medical devices, medical used systems or non medical devices shall be evaluated to the requirements of Clause 16 in the end use application.
- (7) Safety Risk Management: A risk management process complying with the requirements of standard ISO 14971 was inspected during the evaluation of this medical device and/or additional considerations were taken into account since the subject legacy medical device was evaluated to the previous edition of this standard. During the evaluation of this medical device, the risk management decisions affecting the test requirements have been taken into consideration. Changes/Updates in risk management documents that affect the safety of this medical device during its lifecycle shall be communicated to the CSA Group as a condition for continued certification.

MARKINGS

The manufacturer is required to apply the following markings:

- Products shall be marked with the markings specified by the particular product standard.
- Products certified for Canada shall have all Caution and Warning markings in both English and French.

Additional bilingual markings not covered by the product standard(s) may be required by the Authorities Having Jurisdiction. It is the responsibility of the manufacturer to provide and apply these additional markings, where applicable, in accordance with the requirements of those authorities.

The products listed are eligible to bear the CSA Mark shown with adjacent indicators 'C' and 'US' for Canada and US (indicating that products have been manufactured to the requirements of both Canadian and U.S. Standards) or with adjacent indicator 'US' for US only or without either indicator for Canada only.

On the Equipment Exterior:

Equipment is plainly marked in a permanent manner in a place where the details will be readily visible after installation with the following:

- The CSA applicable mark  with optional reference to Standard, CAN/CSA-C22.2 No. 60601-1:14 and ANSI/AAMI ES60601-1:2005/(R)2012, AND C1:2009 AND A2:2010(R)2012
- Manufacturer's identification: Name and/or CSA file number on the same label as the CSA Mark. The name and/or trademark should appear elsewhere on the equipment if only the file number is used on this label.
- Catalogue/Model/Type designation.
- Date of manufacture: Month and year of manufacture or date code. If a serial number is used instead of date of manufacture, a record of serial numbers shall be kept traceable to date of manufacture. (Not related to date of sale).
- Marking on the unit that indicates the manufacturing location if the equipment is manufactured at more than one factory location.
- Complete electrical ratings; in volts (V), hertz (Hz), and amperes (A), Volt-amperes (VA) with the IEC 60417-5032 alternating current symbol adjacent to the marked AC voltage.
- The IEC 60417-5333 "Type BF" symbol for degree of protection against electric shock; only with option "QEVO".

Certificate: 70121971

Project: 80013302

Master Contract: 176164

Date Issued: 2019-10-08

-
- The IEC 417-5021 equi-potential symbol adjacent to the equi-potential earth terminal.
 - The "CAUTION" symbol ISO 7000-0434A on the nameplate, and/or on each or near each output, prefacing CAUTION labels and adjacent to SIP/SOPs.
 - On the power supply cord or on the equipment there is a tag or label indicating that "GROUNDING RELIABILITY CAN ONLY BE ACHIEVED WHEN EQUIPMENT IS CONNECTED TO AN EQUIVALENT RECEPTACLE MARKED 'HOSPITAL ONLY' OR 'HOSPITAL GRADE' " or equivalent wording.

Supplement to Certificate of Compliance

Certificate: 70121971

Master Contract: 176164

*The products listed, including the latest revision described below,
are eligible to be marked in accordance with the referenced Certificate.*

Product Certification History

Project	Date	Description
80013302	2019-10-08	Update cCSAus Certification 70121971/revision 70179363 for Medical Electrical Equipment, Model/Type: KINEVO 900, to cover new software update according to 60601-1 Edition 3.1 without testing based on CPC procedure.
70179363	2018-04-17	Update cCSAus Certification 70121971 for Medical Electrical Component KINEVO 900 to cover Software update from 1.1 to 1.5 according to 60601-1 Edition 3.1.
70154521	2017-09-07	Update of report 70121971 to include Navigation interface and change to Software Version 1.1
70147090	2017-07-23	Update of report 70121971 to include trolley.
70121971	2017-06-27	Initial cCSAus Certification of KINEVO 900.

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
Instruction for Use (IFU)


Dear Madam or Sir,


Please find enclosed the English and Russian Instruction for Use (IFU) for QEVO and QEVO ECU.

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
Address of Record:
Goeschwitzer Strasse 51 - 52
07745 Jena, Germany

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Kurz, Notar

QEVO and QEVO ECU

Digital Exploration Tool and Control Unit - Software Release 1.0

Instructions for Use

Copyright

© 2020, Carl Zeiss Meditec AG, Jena

Registered trademarks

QEVO and QEVO ECU are trademarks or registered trademarks of Carl Zeiss Meditec AG or another company of the ZEISS group in Germany and/or other countries.

All other brands in this document are the property of the respective owner.

Table of Contents

1	Notes on the instructions for use	5
1.1	Product name	5
1.2	Scope.....	5
1.3	Purpose and storage of the documentation	5
1.4	Questions and comments	5
1.5	Conventions in this document	5
1.5.1	Conventions in all text areas.....	6
1.5.2	Conventions in a course of action.....	6
1.6	Other applicable documents	7
2	Safety notes	9
2.1	Target group	9
2.2	Area of use	9
2.2.1	Intended use	9
2.2.2	Indications.....	9
2.2.3	Normal use.....	10
2.3	Responsibilities and duties of the operator	11
2.3.1	Network connection requirements	12
2.3.2	Messages to manufacturer and authorities	12
2.4	Measures and duties of the operator.....	12
2.5	Electromagnetic compatibility	13
3	Product description	15
3.1	Labeling.....	15
3.1.1	Labeling on the product	15
3.1.2	Labeling on the packaging	18
3.1.3	Labeling on the connector panel of the control unit	22
3.2	Design of the product.....	25
3.3	Control unit	26
3.4	Functional description	26
4	Before every use.....	27
4.1	Sterile reprocessing	27
4.2	Visual inspection of the product.....	27
4.3	Configure KINEVO 900 for use of the product.....	27
4.4	Connect the product to the control unit	28
4.5	Function test	30
4.6	Visual inspection of the connector panel of the control unit	31
4.7	Mount the product on the supporting arm system	32

4.8	Incorporate product in navigation system	33
5	Operation	35
5.1	Working with the product	35
5.2	Switching the product off and precleaning it	37
6	Reprocessing	39
7	Maintenance	41
7.1	Maintenance concepts	41
7.1.1	Maintenance via replacement concept	42
7.1.2	Direct maintenance	43
7.2	Implement technical safety test	44
8	Troubleshooting	45
8.1	Localizing faults	45
8.1.1	Response to faults with messages	45
8.1.2	Send the product in for repair	46
8.1.3	Service information	46
8.2	Faults on the product	46
9	Technical specifications	49
9.1	Regulatory Information	49
9.2	Dimensions and weights	49
9.3	Electrical data	49
9.4	Light source	49
9.5	Ambient requirements for operation	50
9.6	Ambient requirements for transport and storage	50
9.7	Guidelines and manufacturer's declaration for electromagnetic compatibility	50
9.7.1	EMC Electromagnetic compatibility IEC 60601-1-2: 2007 (3rd Edition)	50
9.7.2	EMC Electromagnetic compatibility IEC 60601-1-2: 2014 (4th Edition)	55
10	Accessories and spare parts	59
11	Disposal	61
11.1	Disposal of the original packing	61
11.2	Disposal of the product and control unit	61
	Glossary	63
	Index	65

1 Notes on the instructions for use

1.1 Product name

QEVO is a hand-held digital exploration tool which is designated as "product" in these Instructions for Use.

QEVO ECU is the so-called "control unit" for the exploration tool.

1.2 Scope

The present Instructions for Use apply to QEVO and QEVO ECU with Software Release 1.0 and the following identification:

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These Instructions for Use explain the safety features, functions, and performance parameters of the product. They help you to use the product safely and describe proper measures for care and maintenance of the product.

The correct operation of the product is vital for safe and successful operation.

- ▶ Read these Instructions for Use before commissioning and using the product.
- ▶ Keep the Instructions for Use where they are accessible for all users at all times.
- ▶ Give the Instructions for Use to the subsequent owners of the product.

1.4 Questions and comments

- ▶ If you have any questions or information on these Instructions for Use or the software, please contact ZEISS Service.

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This is a cross-reference: Notes on the instructions for use.

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- "File / Save as"
- "My documents / Documents"

Keys, buttons, knobs, levers, and other operating controls are marked by square brackets:

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1.6 Other applicable documents

Document type	Document title
Short Instructions for Use	G-30-1977 QEVO
Reprocessing instructions	G-30-1989 QEVO
Instructions for Use	Valid KINEVO 900 Instructions for Use

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2.2.1 Intended use

QEVO

The product is intended for viewing internal surgical sites during general surgical procedures and for use in visualization of ventricles and structures within the brain during neurological surgical procedures.

QEVO ECU

The product is intended to operate an endoscopic device and to provide a high resolution video signal to a monitor.

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Indication of use	
Discipline:	■ Neurosurgery and spine surgery ■ Orthopedic surgery ■ Oral and maxillofacial microsurgery ■ P&R procedure ■ Not applicable for ophthalmology
Indication:	The QEVO is intended for viewing internal surgical sites during general surgical procedures and for visualizing ventricles and structures within the brain during neurological surgical procedures. The QEVO ECU is intended to operate an endoscopic device and to provide a high resolution video signal to a monitor.

Indication of use	
Contraindication:	The product is contraindicated for treatments in which the endoscopic method is contraindicated.
Patient target group:	No restriction
Age:	No restriction
Gender:	All

2.2.3 Normal use

QEVO

The product has been developed for the use by qualified surgeons to support diagnosis and treatment during general surgical procedures and for use in visualization of ventricles and structures within the brain during neurological surgical procedures. Further fields of application include ENT, spinal and plastic reconstructive surgery.

QEVO ECU

The product exerts its function integrated within a surgical microscope using the power supply, the user interface and the monitor of a surgical microscope.

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Operating personnel

The product may be operated only by instructed and trained persons.

Action

- ▶ Make sure that the operating personnel are appropriately trained and instructed.
- ▶ Make sure that the operating personnel have read and understood the Instructions for Use.
- ▶ Keep the Instructions for Use available at all times for the operating personnel.
- ▶ To simplify access for all operating personnel, order additional copies of the Instructions for Use from ZEISS as required.
- ▶ Specify the competencies for handling the product and state who is authorized to perform what tasks.
- ▶ Determine the reporting obligations for malfunction and damage and make them known. Notification of the manufacturer and authorities.
- ▶ Provide the necessary protective clothing.
- ▶ Regularly check that the legal regulations applicable in your country with regard to accident prevention and work safety are being complied with.

Safety inspections

The safety inspections may only be performed by ZEISS Service.

- ▶ High temperatures and pressures during autoclaving may lead to changes in material which can impair image quality. We recommend performing service after reaching 150 sterilization cycles in order to ensure proper functioning and optimal image quality [44].
- ▶ Regular checking is necessary to constantly ensure proper functioning and a long service life. For example, within the framework of a service and maintenance contract.
- ▶ Comply with the specified time limits.
- ▶ Carry out checks according to the extent specified.

The safety inspections of the product and the control unit should at least comprise the following points:

- Availability of the Instructions for Use
- Visual inspection of the product for damage, as well as legibility of the labels
- Functional test of the light, the camera, the image quality and the connector on the product.
- Leakage current test

- Test of protective ground conductor
- Functional test of the socket for the product on the connector panel of the control unit.

Maintenance and Inspection

- ▶ In order to constantly ensure proper functioning of the device and reach its expected useful life: Comply with the maintenance and inspection intervals specified in these Instructions for Use.

Modifications to the product

Modification of the product is not allowed. This product must not be modified without the manufacturer's approval. If the product is modified, suitable inspections and tests must be performed to ensure that the product can still be used safely.

2.3.1 Network connection requirements

The control unit for the product is permanently installed in the KINEVO 900. The requirements for connecting the KINEVO 900 to a network are applicable here.

Action

- ▶ Please read and follow the Instructions for Use for the KINEVO 900 [7].

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Electrical safety

Action

- ▶ Do not touch the product if your body is electrostatically charged and the product is not grounded.
- ▶ Please note the information on electromagnetic compatibility (EMC).
- ▶ Live parts are freely accessible in the interior of the product and the control unit. If you remove the housing, you run the risk of electric shock. Never open the housing!

The electrical interconnection conditions demand that the applied parts from other ME devices used within the configuration for endoscopic application must be applied parts of the type BF or CF.

Ambient conditions

- ▶ Operate the product only within the limits of the prescribed ambient conditions [50].
- ▶ Ensure that the installation conditions and operation meet surgical requirements:
 - Low vibration
 - Clean environment
 - Avoid extreme mechanical stress
- ▶ Do not use power-operated products included in the delivery package
 - in explosive atmospheres,
 - at a distance of less than 25 cm from flammable anesthetics or volatile solvents such as alcohol, benzine or similar substances.

Symbols and labels

- ▶ Note the symbols and labels attached to the product!

Transport

- ▶ Transport the device over long distances (e.g. relocation, return for repair) only in its original packaging or special return packaging.
- ▶ Please contact your dealer or ZEISS Service for this purpose.

2.5 Electromagnetic compatibility

The product is subject to specific requirements with regard to electromagnetic compatibility (EMC). The following factors can cause EMC malfunctions:

When the product is used, observe the following EMC precautions regarding the EMC:

- Portable and mobile RF communication devices in the vicinity of the product.
- Other products which are set up in the vicinity.
- Accessories and cables and spare parts that are not specified in these Instructions for Use and not sold by ZEISS as spare parts.

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- ▶ Observe the Instructions for Use.
- ▶ Observe the EMC guidelines in the section "Specifications".
- ▶ If you use radios or wireless transmission components: Maintain a minimum distance of 30 cm to all components of the product.
- ▶ If you set up the product in the vicinity of other devices: Check whether normal operation is possible in this arrangement.

Procedure

Action

Empty page, for your notes

3 Product description

3.1 Labeling

CAUTION!

Function failure due to incorrect operation and wrong storage!

If the labeling on the product and connector panel of the control unit is not correctly understood or read, the product can be operated and stored incorrectly.

- ▶ Make sure that all labeling on the product and connector panel of the control unit is legible.
- ▶ When operating the product, observe all labeling on the product and connector panel of the control unit.
- ▶ Observe all labeling on the packaging for storage and transportation.

3.1.1 Labeling on the product

3.1.1.1 Markings on the Upper Side

Figure 1: Product

Pos.	Symbol	Name	Explanation
1	45°	Degree specification	Direction of view

3.1.1.2 Labeling on the lower surface



Figure 2: Product

Pos.	Symbol	Name	Explanation
1	INT 500s/1000s	Intermittent mode	Operating mode: Do not use [26] the product in continuous operation.
2		UDI	UDI labeling
3		Part number	Reference number
4		CE mark	-
5		Applied part BF	Fulfills the BF conditions for protection against electric shock according to IEC 60601-1.
6		Serial number	-
7		Maintenance interval	Informs about implemented maintenance intervals.

3.1.1.3 Labeling on the connector

Figure 3: Connector

Pos.	Symbol	Name	Explanation
1		Arrow	Insertion direction on the connector panel of the control unit QEVO ECU
2		Applied part BF	Satisfies the BF conditions in accordance with IEC 60601-1 for protection against electrical shock.

3.1.2 Labeling on the packaging

3.1.2.1 Packaging sign






Figure 4: Packaging label

Pos.	Symbol	Name	Explanation
1		Keep away from sunlight	Protect the packaging and packaged contents from heat and solar radiation.
2		Keep dry	Protect packaging and packaged contents from moisture.
3		Instructions for Use	Read and observe the Instructions for Use.
4		Non-sterile	The product in the packaging has been cleaned, but it is not sterile. The products must be completely conditioned before use.
5		-	The product may only be used by a physician or by another person on behalf of a physician.

Pos.	Symbol	Name	Explanation
6		Unsuitable for MR	The product may not be operated or stored in an environment in which magnetic resonance occurs.
7		UDI	UDI labeling
8		UDI Identity No.	Contains information on the unique identification of the product.
9		CE mark	-
10		Endoscope	Type of product
11	Made in Germany	Country of manufacture	The product was manufactured in Germany.
12		Packing unit	Packaging contains one product.
13		Air pressure limit	The product may only be transported and stored at a minimum atmospheric pressure of 500 hPa and a maximum atmospheric pressure of 1060 hPa.
14		Humidity limit	The product may only be transported and stored at a minimum humidity level of 10% RH and a maximum humidity level of 92% RH.
15		Temperature limit	The product may only be transported and stored at temperatures no lower than -20 °C and no higher than +60 °C.
16		Serial number	-
17	QEVO	Product name	-
18		Part number	Reference number


Pos.	Symbol	Name	Explanation
19		Manufacturing information	Contains the name and address of the manufacturer and the manufacture date of the product.
20	INT 500s/1000s	Intermittent mode	Operating mode: Do not use [26] the product in continuous operation.

3.1.2.2 Packaging

Figure 5: Packaging

Pos.	Symbol	Name	Explanation
1		Indication of direction for transporting and storing the packaging "top"	Arrows must point upwards when placing the package.
2		Fragile, handle with care	The packaged contents are fragile.
3		Keep dry	Protect the packaging and packaged contents from moisture.

Pos.	Symbol	Name	Explanation
4		Humidity limitation	The product may only be transported and stored at an air humidity of min. 10 % and max. 92 % RH.
5		Atmospheric pressure limitation	The product may only be transported and stored at an atmospheric pressure of min. 500 hPa and max. 1060 hPa.
6		Temperature limit	The product may only be transported and stored at a temperature range of min. -20 °C and max. +60 °C.
7		Number of packaging units	Packaging contains one product.

3.1.3 Labeling on the connector panel of the control unit

Connector panel

The connector panel of the control unit for the product is located on the KINEVO 900.



Figure 6: Connector panel of the control unit

1		Endoscope	Type of product
2		Rating label	Contains the performance data of the control unit [23]
3		Product label	Contains specific information about the product [24]

3.1.3.1 Power rating plate of the control unit


Figure 7: Rating plate of the control unit

Pos.	Name	Explanation
1	Carl Zeiss Meditec AG	Manufacturer
2	Goeschwitzer Strasse 51 - 52 07745 Jena, Germany	Manufacturer's address
3	Max. 48VA	Connected load
4	CE mark	-
5	CSA label	This medical device satisfies the US American and Canadian safety regulations for medical electrical systems.
6	Observe disposal regulations	Do not dispose of electrical or electronic devices along with normal domestic waste. For more information on the disposal of electrical and electronic devices, see the section entitled "Disposal".
7	MADE IN GERMANY	Country of manufacture
8	24 V DC	Rated voltage
9	QEVO ECU	Product name of control unit

3.1.3.2 Product identification label


Figure 8: Product label

Pos.	Name	Explanation
1	MD	Identifies the product as a medical device
2	UDI	UDI labeling
3	UDI Identity No.	Contains information about the manufacturing date and the consecutive serial number
4	Date of manufacture	Year-Month-Day
5	SN	Serial Number
6	REF	Catalog number

3.2 Design of the product


Figure 9: Product overview

1	Shaft*	2	Hand grip**
3	Cable	4	Connector

* The shaft contains rod lens optics: light is conducted through rod lenses made of quartz glass and broken by air lenses between the rods.

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** The hand grip contains a camera and a light source.

3.3 Control unit


Figure 10: Control unit

1	Control unit	2	Connection socket for product connector
3	Connector panel	4	Ventilation grille for control unit

The control unit for the product is permanently installed in the KINEVO 900. Only the connector panel for the connector of the product is visible.

3.4 Functional description

If the product is connected to the connection socket of the control unit, it is ready for operation within a few seconds. The light source in the hand grip illuminates the surgical field. The camera integrated in the product takes pictures or videos and displays them on the monitor.

The electronics heat the handle. The following measures prevent overheating of the handle:

- Put the product aside. When the automatic light switch-off function is activated, the light source automatically switches off after approx. 40 seconds, allowing the handle to cool down.
- 500s/1000s intermittent mode: Operate the product for a maximum of 500 seconds (eight minutes) and then switch the product off for at least 1000 seconds (17 minutes).
- The control unit automatically decreases the light intensity and camera performance of the product as its temperature increases.

4 Before every use

4.1 Sterile reprocessing

Action

1. The product must be reprocessed in sterile condition before each use, including before it is used for the first time.
2. Please read and observe the Reprocessing Instructions [7].

4.2 Visual inspection of the product

CAUTION!

Functional failure during operation!

The product must not be used if it shows signs of damage.

- Perform a visual inspection prior to each use.
- Keep a spare product available for use in case the product is damaged.

Prerequisite

- ☒ The product has been reprocessed in sterile condition
- ☒ Sterile gloves
- ☒ Sufficient ambient light

Action

1. Wear sterile gloves to avoid any contamination of the product.
 2. Inspect the product: Its shaft, hand grip, cable and connector must not have any deposits, cracks or defects.
 3. Check the tip of the shaft: The tip of the shaft must not have any cracks or deposits and must not be fogged. The glass surface must be clean and smooth.
 - ⇒ Do not use the product if there is any visible contamination. Do not use the product if it has sharp corners or edges, protrusions or rough or scratched surfaces that could injure the patient.
 4. Perform a visual inspection before each use of the product.
- ZEISS recommends keeping a spare product available for use in case the product is not fit for operation.

4.3 Configure KINEVO 900 for use of the product.

In order to operate the product via the KINEVO 900, several adjustments must first be performed on the KINEVO 900.

These adjustments concern the hand grip, the foot control panel and the software of the KINEVO 900 and are not components of the product.

Action

- Please read and follow the Instructions for Use for the KINEVO 900 [7].

4.4 Connect the product to the control unit

A sterile and a non-sterile person are required to connect the product. The sterile person places the sterile product near the patient. The non-sterile person plugs the product cable into the connection socket on the control unit.

WARNING!

Risk of explosion!

There is an increased risk of fire or explosion in oxygen-enriched environments.

- ▶ Do not operate the product or the control unit in explosive atmospheres or at a distance of less than 25 cm from flammable anesthetics or volatile solvents such as alcohol, benzine or similar substances.

CAUTION!

Tripping hazard!

The cable between the product and the control unit in the KINEVO 900 increases the risk of tripping.

- ▶ Lay the cable so that it does not obstruct the workflow in the operating room.

CAUTION!

Connection of inadmissible or defective accessories will result in an increased leakage current!

If you plug unapproved or faulty accessories into the connection socket on the control unit, this may result in increased voltages which are then transmitted to the patient.

- ▶ When setting up your system, observe the (German) Medical Device Operator Ordinance or other applicable national legislation.

NOTE

No camera image is displayed on the monitor!

If you connect a product from another manufacturer to the connection socket of the control unit, no camera image will be displayed.

- ▶ Only plug compatible products into the connection socket on the control unit.
- ▶ Read and observe the Instructions for Use pertaining to these products.

Prerequisite

- ☑ KINEVO 900 has been cleaned and disinfected, see Instructions for Use for the KINEVO 900 [7]
- ☑ The product has been activated on the KINEVO 900, see Instructions for Use for the KINEVO 900 [7].
- ☑ The product has been reprocessed in sterile condition
- ☑ The product is lying in the sterilized Instru-Safe instrument tray.

Action

1. Sterile person: Take the product with the cable out of the Instru-Safe instrument tray.
2. Sterile person: Place the product safely in the sterile area.
3. Unsterile person: Plug the connector into the connection socket of the control unit.

⇒ The product is ready for operation immediately.

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4.5 Function test

CAUTION!

Functional failure during operation!

The product must not be used unless it is fully functional.

- ▶ Perform a functional test prior to each use.

CAUTION!

Misdiagnosis due to monitor display and video images!

The monitor on the KINEVO 900 and the camera in the product are neither calibrated nor designed for diagnostic purposes. Reproduced photos and videos may contain deviations regarding their scale, shape, color and contrast.

- ▶ Observe a structure in the surgical field through the eyepiece of the KINEVO 900 and compare this view with the camera image of the product. The camera image must have a similar quality.
- ▶ Do not use the displayed videos and photos for diagnostic purposes.

CAUTION!

Risk of glare!

Looking directly into the light source at the shaft tip may result in damage to the eye.

- ▶ Do not look into the shaft tip once the product is connected to the connection socket on the control unit and ready for operation.

NOTE

Functional failure during operation!

The product must not be used unless it is fully functional.

- ▶ Keep a spare product available for use in case the product is not fit for operation.

Prerequisite

- ☒ Check that the KINEVO 900 is switched on and ready for use.
- ☒ The product is connected to the connection socket of the control unit.
- ☒ The product is activated on the KINEVO 900, see Instructions for Use for the KINEVO 900. [7]
- ☒ A button on the handgrip and/or on the foot control panel of the KINEVO 900 has been configured for switching the product on, see Instructions for Use for the KINEVO 900. [7]

Action

1. On the monitor of the KINEVO 900, activate the product's user interface.
2. Press the [QEVO] button on the hand grip or the foot control panel.
 - ⇒ The product's light and camera switch themselves on.
 - ⇒ The product's camera image is displayed as a full-screen image. The camera image of KINEVO 900 is displayed on the touchscreen as Picture in Picture.
3. Make sure that all functions can be adjusted on the product.
4. Point the tip of the shaft at a sterile object located nearby.
5. Make sure that you see a sharp, bright display in good image quality.
 - ⇒ Do not use the product if you perceive streaking, shifts in color or image flickering or cannot produce a sharp, bright picture in good image quality.
6. Perform this functional test prior to each use of the product.

4.6 Visual inspection of the connector panel of the control unit

Fans cool the control unit during operation.

NOTE

Functional failure during surgery!

If the vent slots are covered, the control unit overheats and fails to operate.

- Make sure before and during every use that the vent slots on the connector panel of the control unit are not covered.

Action

1. Make sure that the vent slots on the connector panel of the control unit are not covered.
2. Conduct a visual inspection of the connector panel before every use of the product.

4.7 Mount the product on the supporting arm system

The operator is responsible for using the product with a supporting arm system. The operator must ensure that the manufacturer of the supporting arm system describes it as suitable for the task.

WARNING!

Injury to tissue!

If the operating table or the product is moved inadvertently, the shaft can touch the patient and cause tissue trauma.

- ▶ Use the product in connection with a suitable supporting arm system only if the patient and the operating table are fixed in place so that nothing can move, and observe all the manufacturer's data.
- ▶ Use an appropriate clamping device that tightly holds the product in place so that the product cannot slip, and observe all the manufacturer's data.

CAUTION!

Risk of injury due to falling product!

If an unsuitable supporting arm system and clamping devices are used for the product, the product can fall into the surgical field and injure the patient.

- ▶ Use only a suitable supporting arm system with a clamping device suitable for the product, and observe all the manufacturer's data.

Prerequisite

- ☒ The suitable supporting arm system is properly assembled; see the manufacturer's Instructions for Use.
- ☒ A clamping device suitable for the product is available; see the manufacturer's Instructions for Use.
- ☒

Action

1. Mount the product on the supporting arm system with the clamping device.
2. Position the product on the patient.
3. Switch off the product's automatic light switch-off.

For more information on the supporting arm system, refer to the manufacturer's instructions for use.

4.8 Incorporate product in navigation system

The operator and the manufacturer of the navigation system are responsible for use of the product with a navigation system.

CAUTION!

Incorrect reference in the navigation system!

Surgical fields will be wrongly calculated if the references are incorrect.

- ▶ Only allow properly instructed and trained personnel to incorporate the product into the navigation system.
- ▶ Read and observe the manufacturer's data on the navigation system.
- ▶ Use an appropriate clamp that firmly encloses the product for the navigation probe, so that the navigation probe cannot slip. Observe the manufacturer's data.

Prerequisite

- ☒ A clamp that is suitable for the product is available.

Action

- ▶ Mount the navigation probe on the product securely with a suitable clamp.

For more information on the navigation system, refer to the manufacturer's Instructions for Use.

Empty page, for your notes

5 Operation

More detailed information on the product's user interface and settings is provided in the Instructions for Use of the KINEVO 900.

Action

- Read and follow the Instructions for Use for the KINEVO 900 [7] before proceeding to work with the product.

5.1 Working with the product

WARNING!

The use of a defibrillator can cause injury to tissue!

If complications arise during the application that require the use of a defibrillator and the shaft of the product is in the surgical field, the shaft touches the tissue and a tissue trauma develops.

- If complications arise and a defibrillator is used, you must remove the shaft from the surgical field.

WARNING!

Injury to tissue due to contact with shaft!

Only tissue located in the camera's field of view is visible on the monitor, i.e. the tissue surrounding the entire surgical field is not displayed. If the shaft touches the surrounding tissue, this causes tissue trauma.

- Guide the shaft carefully and cautiously within the surgical field.
- If you move the product within the surgical field, do not look at the camera image; instead, move the product under direct visual control.
- If you move the product within the surgical field, avoid touching any tissue.
- Do not operate the product near a magnetic field (e.g. an MR scanner), as magnetic fields can cause unintentional movement of the product.

CAUTION!

Electric shock due to increased leakage current!

A voltage being able to injure patient and user is applied on plug connectors.

- Never touch the contacts of plug connectors during contact with the patient.

⚠ CAUTION!

Misdiagnosis due to monitor display and video images!

The monitor on the KINEVO 900 and the camera in the product are neither calibrated nor designed for diagnostic purposes. Reproduced photos and videos may contain deviations regarding their scale, shape, color and contrast.

- ▶ Do not use the displayed videos and photos for diagnostic purposes.

⚠ CAUTION!

Injury of tissue due to unsupervised light source!

If the product is used with a high light intensity over a longer period of time, tissue damage can occur in the area of the illuminated surgical field.

- ▶ Switch the automatic switch-off function on.
- ▶ Do not leave the product unattended while it is switched on.
- ▶ Operate the product in intermittent mode.
- ▶ Do not put the product down near the patient.

⚠ CAUTION!

Heat build-up on the handle!

The handle heats up during continuous operation.

- ▶ Do not operate the product in continuous operation, but rather in intermittent mode.

⚠ CAUTION!

Risk of glare!

Looking directly into the light source at the shaft tip may result in damage to the eye.

- ▶ Do not look into the shaft tip when the product is connected to the connection socket of the control unit and ready for operation.

NOTE

Failure of light on the product!

The light emitted from the tip of the shaft may fail due to a system crash.

- ▶ Hold the product steady.
- ▶ Use the light on the KINEVO 900 or your surgical lamps to move the shaft out of the surgical field, if possible without touching the surrounding tissue.

NOTE

Functional failure due to harmful environmental conditions!

The product must not be used in the event of a functional failure.

- ▶ Operate the product under the specified environmental conditions [50] only.

Prerequisite

- ☑ The functional test has been completed [30].
- ☑ KINEVO 900 has been switched on and is ready for operation.
- ☑ The product has been activated on the KINEVO 900, see Instructions for Use for the KINEVO 900 [7].
- ☑ A button on the handle and/or on the foot control panel of the KINEVO 900 has been configured for switching the product on, see Instructions for Use for the KINEVO 900 [7].
- ☑ White balance has been performed, see Instructions for Use for the KINEVO 900 [7].

Action

1. Pick up the product.
2. Press a [QEVO] function key on the foot control panel or the handle of the KINEVO 900 to activate the product.
3. Slowly and carefully guide the shaft into the surgical field. Guide the product under visual control so that you can see exactly where the tip of the shaft is located. To the extent possible, make sure that the shaft does not touch any tissue.
4. If necessary, adapt the orientation of the product's camera image to the viewing direction of the KINEVO 900.
5. You can put the product down on a sterile surface at any time during the operation.
 - ⇒ If the automatic switch-off function is activated, the product automatically switches off after approx. 40 seconds.
6. If the displayed image is too dark or too bright, use the yellow function buttons on the right handle of the KINEVO 900 to adjust the brightness.
7. Then operate the product in intermittent mode: Use the product for a maximum of eight minutes and then switch the product off for at least 17 minutes to prevent the handle from overheating.

5.2 Switching the product off and precleaning it

CAUTION!

Heat build-up on the handle!

The handle heats up during operation.

- Let the product cool down briefly before cleaning it.

NOTE

Functional failure during operation!

The product must not be used if it shows signs of damage.


- Use the supplied Instru-Safe instrument tray to avoid damage to the product during processing.
Read and observe the attached reprocessing instructions for the "QEVO" (G-30-1989) [7].

Prerequisite

- ☒ Enzymatic cleaning solution, e.g. Cidezyme (ASP®), prepared according to the manufacturer's Instructions for Use
- ☒ Two dry, lint-free cloths
- ☒ Instru-Safe instrument tray
- ☒ Sterile gloves

Action

1. Switch the product off by pressing the [QEVO] button on the handle or the foot control panel of the KINEVO 900, see Instructions for Use for the KINEVO 900 [7].
2. Disconnect the connector of the product from the connection socket of the control unit.
3. Let the product cool down briefly before cleaning it.
4. Moisten a cloth with enzymatic cleaning solution.
5. Wipe off the product, the cable and the connector with the moistened but not dripping wet cloth until any visible surgical residues are removed as completely as possible, but at least for one minute.
6. Wipe off the product, the cable and the connector with a dry, lint-free cloth.
7. Place the product with cable and connector in the Instru-Safe instrument tray. Read and observe the attached "QEVO" Reprocessing Instructions (G-30-1989).

8. Send the product in for reprocessing and make sure that it is reprocessed within six hours.

6 Reprocessing

QEVO must be replaced if the image is clouded, distorted, or there is no image at all.

QEVO can - with appropriate care and, if undamaged and not contaminated – be used up to 150 times. Any further reuse or the use of damaged and / or contaminated products is the responsibility of the user.

WARNING!

Contamination with pathogenic microorganisms!

If the product is reprocessed using unsuitable methods or cleaning agents, pathogenic microorganisms can settle on the product.

- ▶ Reprocess the product only according to the validated methods for this product.
- ▶ Only allow briefed and trained personnel to reprocess the product.
- ▶ Only use compatible cleaning agents when reprocessing the product.
- ▶ Dispose of the product as soon as there is any suspicion of prion contamination. Observe the respective prevailing national regulations for the disposal.

Action

1. Reprocess the product according to the validated processes described in the attached reprocessing instructions "QEVO" (G-30-1989).
2. Please be aware, that the number of plug cycles for a QEVO can differ from the number of sterilization cycles. The tracking of the sterilization cycles is in the responsibility of the user. The counting of plug cycles is a reference indicator for your convenience.

Empty page, for your notes

7 Maintenance

The product is maintained exclusively by your local ZEISS Service. The control unit is permanently installed in the KINEVO 900, and undergoes a regular safety inspection at the same time as the KINEVO 900.

CAUTION!

Increased leakage current!

Improperly performed maintenance and service work on the product and control unit can disconnect lines and plug connections. Housing parts are live.

- You should only have the product and control unit serviced by ZEISS Service.

NOTE

Functional failure!

Improperly performed maintenance and service work on the product and control unit can lead to functional failure.

- You should only have the product and control unit serviced by ZEISS Service.

7.1 Maintenance concepts

QEVO must be replaced if the image is clouded, distorted, or there is no image at all.

QEVO can - with appropriate care and, if undamaged and not contaminated – be used up to 150 times. Any further reuse or the use of damaged and / or contaminated products is the responsibility of the user.

NOTE

Ensuring the functioning and image quality through maintenance:

We recommend performing service after reaching 150 sterilization cycles in order to ensure proper functioning and optimal image quality.

- Comply with the specified maintenance intervals.

You have two options for having your product maintained by ZEISS Service:

1. Replacement concept

- You receive a replacement product and send in [42] your own product for maintenance.
- The replacement product has a different serial number.
- The replacement product is marked in the "Maintenance interval" [16] field.

2. Direct maintenance

- You send your product in for maintenance and it is returned to you again [43] completely repaired and reprocessed.
- Maintenance can be performed on your product twice. The product must be replaced after it is undergone maintenance twice.

7.1.1 Maintenance via replacement concept

Maintenance via the replacement concept means that you get a completely repaired and reprocessed product and only then need to send in your own product for maintenance. Your advantage here is that you are provided with a replacement product within a short period of time.

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If products which have not been reprocessed are returned to ZEISS Service, pathogenic germs may spread. Cross-contamination can occur.

- Reprocess the product completely before sending it in to ZEISS Service for maintenance.

Prerequisite

- ☒ The QEVO image is cloudy, unsharp or noisy or no image at all can be seen

We recommend performing service after reaching 150 sterilization cycles in order to ensure proper functioning and optimal image quality.

Action

1. Get in touch with your ZEISS Service organization and inform them about the maintenance of your product according to the "replacement concept".
 - ⇒ You will receive a replacement product. The original packaging of the replacement product also includes a new bag, a certificate of hygienic compliance, a certificate for repair/service and reprocessing and a return shipment seal.
2. Remove the replacement product from the original packaging.
3. Completely reprocess the product that must be sent in for maintenance.
4. Place the reprocessed product in the supplied bag and seal it.
5. Place the product packaged in the bag in the original packaging.
6. Fill out the certificate of hygienic compliance and place it in the original packaging.
7. Seal the original packaging and affix the supplied seal to it.
8. Return the package to your ZEISS Service organization.

You can find the ZEISS contact partner for your country on the Internet in the following website: www.zeiss.com/med.

7.1.2 Direct maintenance

Direct maintenance means that you send your product in for maintenance and the same product is returned to you completely repaired and reprocessed. This maintenance process takes longer than maintenance via the replacement concept.

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If products which have not been reprocessed are returned to ZEISS Service, pathogenic germs may spread.
Cross-contamination can occur.

- Reprocess the product completely before sending it in to ZEISS Service for maintenance.

Prerequisite

- ☒ The QEVO image is cloudy, unsharp or noisy or no image at all can be seen

We recommend performing service after reaching 150 sterilization cycles in order to ensure proper functioning and optimal image quality.

Action

1. Get in touch with your ZEISS Service organization and inform ZEISS Service about the maintenance of your product according to the "direct maintenance" concept.
2. If you no longer have the original packaging with an unused bag, seal, certificate of hygienic compliance and certificate for repair/service and reprocessing, order a new packaging unit [59].
3. Reprocess your product completely.
4. Remove the unused bag from the original packaging.
5. Place the reprocessed product in the bag and seal it.
6. Place the product packaged in the bag in the original packaging.
7. Fill out the certificate of hygienic compliance and place it in the original packaging.
8. Seal the original packaging and affix the seal to it.
9. Return the package to your ZEISS Service organization.

You can find the ZEISS contact partner for your country on the Internet in the following website: www.zeiss.com/med.

7.2 Implement technical safety test

The safety inspection serves to determine and evaluate product safety. The user of this product is obligated to have safety inspections performed and documented in compliance with IEC 62353.

CAUTION!

Risk of injury without safety inspection!

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- ▶ Have a safety inspection in compliance with IEC 62353 performed on time and to the prescribed extent. At the same time, be sure to observe the relevant national regulations.
Interval for the product: After 150 sterilization cycles or, at the latest, every two years.
Interval for the control unit: Based on the interval for the KINEVO 900.

Prerequisite

- ☒ Only persons certified by ZEISS may perform the safety inspection.

Action

- ▶ Send the product in [41] to ZEISS Service after 150 sterilization cycles or, at the latest, every two years and have a safety inspection performed and documented.
- ▶ Have a safety inspection performed and documented on the connection socket and on the control unit according to country-specific requirements and including a leakage current test concurrent to the safety inspection on the KINEVO 900. You will find detailed information on the safety inspection in the Instructions for Use of the KINEVO 900.

8 Troubleshooting

What to do in the event of malfunctions

CAUTION!

Increased leakage current!

Incompetently conducted maintenance and service work on the product and control unit can detach lines and plug connectors. Housing parts are under voltage.

- ▶ Have the product and control unit serviced only by ZEISS Service.

NOTE

Functional failure!

Incompetently conducted maintenance and service work on the product and control unit can lead to functional failure.

- ▶ If an error occurs that you cannot rectify on the basis of this section, then label the product as non-functioning.
- ▶ Inform ZEISS Service.

8.1 Localizing faults

If a particular event occurs during operation or the control unit registers an irregularity, the control unit indicates this fault by displaying a message on the monitor.

Faults are saved in a service file. You can export this service file and transmit it to ZEISS.

8.1.1 Response to faults with messages

Action

1. If an error message appears on the monitor, read through the message text carefully.
2. Tap on the error message to close it.
3. Correct the displayed error/fault.
4. If this error message remains displayed on the monitor, export the service file of the KINEVO 900 and send it to ZEISS Service. Read and observe the Instructions for Use of the KINEVO 900.
5. If the error messages prompt you to replace the product, send the defective product to ZEISS Service [46] for repair.

8.1.2 Send the product in for repair

Send your product to ZEISS Service for repair if

- a fault occurs which you cannot remedy based on the section Product faults [46] or
- error messages appear on the monitor which prompt you to replace the product.

ZEISS offers you two options for having your product repaired by ZEISS Service:

- **Replacement concept:** You receive a fully repaired and reprocessed replacement product and only then send your own product in for repair.
- **Direct maintenance:** You send in your product for repair and receive the same product back again, completely repaired and reprocessed, after a longer period of time.

The procedure for sending in the product for repair is based on the same principle as sending in the product for maintenance:

Action

- ▶ Inform ZEISS Service about the repair of your product according to the "Replacement concept [42]".
- ▶ Inform ZEISS Service about the repair of your product according to the "Direct maintenance [43]" concept.

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8.2 Faults on the product

Fault	Cause	Remedy
Image is cloudy	Glass surface in shaft tip is dirty.	▶ Clean glass surface in shaft tip. ▶ Check water quality. ▶ Then reprocess the product completely.
	Lens system in shaft tip is leaky or defective.	▶ Send product in [46] to ZEISS Service for repair.
Image is impaired or drops out completely	Radio communications such as DECT, WLAN, Bluetooth, cell phones, RFID transmitters and other emitting electronic devices can disturb the product.	▶ Read and observe the guidelines and the manufacturer's declaration concerning electro-magnetic compatibility [50]. ▶ Remove all electronic devices located in the vicinity of the product. ▶ Restart the product by unplugging its connector from the connection socket on the control unit and then plugging it back in.

Fault	Cause	Remedy
The image is too dark, the illumination is too low	Image brightness set too low	► Increase the image brightness., see Instructions for Use for the KINEVO 900 [7].
	Glass surface in shaft tip is dirty.	► Clean glass surface in shaft tip. ► Check the quality of your cleaning water. ► Then reprocess the product completely.
	Light source is defective.	► Send product in [46] to ZEISS Service for repair.
Image is color-cast.	White balance was not performed.	► Perform white balance, see Instructions for Use for the KINEVO 900 [7].
	Glass surface in shaft tip is dirty.	► Clean glass surface in shaft tip. ► Check water quality. ► Then reprocess the product completely.
	Light source is defective.	► Send product in [46] to ZEISS Service for repair.
Test image appears.	QEVO is inserted incorrectly or not at all.	► Check the plug connection; replug the connector into the connection socket of the control unit if necessary.
Surfaces of the product are corroding, exhibit staining or discolorations.	The product was insufficiently cleaned (e.g. remaining proteins).	► Rub the product down with a cleaning solution and dry it thoroughly with a lint-free cloth. ► Clean the product manually. ► Then reprocess the product completely.
	The product was rinsed insufficiently between the reprocessing phases, especially prior to sterilization.	► Rinse the product sufficiently between the reprocessing phases.
	The chloride concentration of the water is too high.	► Check water quality.
	The concentration of minerals (e.g. lime) or organic substances in the water is too high.	► Check the quality of your cleaning water. ► Use only demineralized water if necessary.
Surfaces of the product are corroding, exhibit staining or discolorations.	Water or sterilizing steam contains heavy metal ions and/ or silicates; its iron, copper and manganese content is increased.	► Check the quality of your cleaning water. ► Use only demineralized water if necessary.
	The cleaning and disinfectant solution is contaminated.	► Renew the cleaning and disinfectant solution regularly.
	Extraneous rust, e.g. due to steam containing rust or reprocessing together with damaged or rusty instruments	► Check the supply systems. ► In the case of simultaneous reprocessing of several products, check for material compatibility and prior damage. ► Avoid contact of the product with other products.

Empty page, for your notes

9 Technical specifications

9.1 Regulatory Information

Classification of the Product and the Control Unit according to IEC 60601-1

QEVO and QEVO ECU are classified as follows:

- QEVO ECU fulfills the level of protection against electric shock: Class 1
- QEVO fulfills the type BF conditions for protection against electric shock
- Mode of operation QEVO: Intermittent
- Mode of operation QEVO ECU: Continuous

Classification of the Product and the Control Unit according to IEC 60601-1-2

QEVO and QEVO ECU are classified as follows:

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9.2 Dimensions and weights

Designation	Unit	Value
Weight of hand grip with shaft	g	250
Length of shaft	mm	120
Diameter of shaft	mm	3.6
Length of cable	mm	5000

9.3 Electrical data

Designation	Unit	Value
Operating voltage	VDC	24 ($\pm 10\%$)
Current consumption	A	2

9.4 Light source

Designation	Unit	Value
Luminous flux	LM	20....35

9.5 Ambient requirements for operation

Designation	Unit	Value
Temperature	°C	+10 ... +40
Relative humidity	%	30 ... 75
Atmospheric pressure	hPa	700 ... 1060

9.6 Ambient requirements for transport and storage

Designation	Unit	Value
Temperature	°C	-20 ... +60
Relative humidity	%	10 ... 92
Atmospheric pressure	hPa	500 ... 1060

9.7 Guidelines and manufacturer's declaration for electromagnetic compatibility

9.7.1 EMC Electromagnetic compatibility IEC 60601-1-2: 2007 (3rd Edition)

The product and the control unit are subject to specific precautions with regard to electromagnetic compatibility (EMC). In order to avoid the occurrence of EMC interference, the product may only be installed, operated and maintained in the manner indicated in these Instructions for Use and only with components supplied by ZEISS.

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- ▶ Only use accessories, transformers, cables and spare parts which are specified in these Instructions for Use or which are approved by ZEISS for this device.
- ▶ The product can be disturbed by other devices, also when these other devices conform with the emission requirements valid for them according to CISPR.

Do not operate the product in the direct vicinity of another device. If this cannot be avoided, you must monitor the proper function of the product and the other device.

- Portable and mobile radio communication equipment should not be used at a shorter distance from the product, including cables, than the recommended safety distance that is calculated using the equation applicable to the transmission frequency involved.
- Please note the EMC guidelines on the following pages.

9.7.1.1 Electromagnetic interference

QEVO and QEVO ECU are intended for operation in an electromagnetic environment as specified below. The customer or the user of QEVO and QEVO ECU is responsible for ensuring that both products are operated in such an environment.		
Interference measurements	Compliance	Electromagnetic environment - guidelines
RF emissions as per CISPR 11	Group 1	QEVO and QEVO ECU use RF energy only for their internal functions. Their RF emission is very low and it is improbable that it will interfere with neighboring electronic devices.
RF emissions as per CISPR 11	Class A	The QEVO and QEVO ECU are suitable for use in other establishments than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Not applicable	
Emission of voltage fluctuations/flicker as per IEC 61000-3-3	Not applicable	

9.7.1.2 Electromagnetic immunity for ME equipment and ME systems

QEVO and QEVO ECU are intended for operation in an electromagnetic environment as specified below. The customer or the user of QEVO and QEVO ECU is responsible for ensuring that both products are operated in such an environment.			
Electromagnetic immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment - guidelines
Electrostatic discharge (ESD) as per IEC 61000-4-2	±6 kV contact discharge ±8 kV air discharge	±6 kV contact discharge ±8 kV air discharge	Floors should be made of wood or concrete or be covered with ceramic tiles. If the flooring contains synthetic materials, the relative humidity must be at least 30%.

Fast transient/burst immunity as per IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	The quality of the supply voltage should be that of a typical business or hospital environment.
Surges as per IEC 61000-4-5	± 1 kV line-to-line voltage ± 2 kV line-to-ground voltage	± 1 kV line-to-line voltage ± 2 kV line-to-ground voltage	The quality of the supply voltage should be that of a typical business or hospital environment.
Voltage dips, short interruptions and voltage variations as per IEC 61000-4-11	$< 5\% U_T$ ($> 95\%$ dip in U_T) for 1/2 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $< 5\% U_T$ (95% dip in U_T) for 5 s	$< 5\% U_T$ ($> 95\%$ dip in U_T) for 1/2 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $< 5\% U_T$ (95% dip in U_T) for 5 s	The quality of the supply voltage should be that of a typical business or hospital environment. If the user of QEVO and QEVO ECU require continued function even in the event of interruptions in the power supply, we recommend to power the QEVO and QEVO ECU from an uninterruptible power supply or a battery.
Magnetic field for supply frequency (50/60 Hz) as per IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields in the supply frequency should correspond to the typical values that are found in business and hospital environments.
Note: U_T is the AC supply voltage prior to application of the test level.			

9.7.1.3 Electromagnetic immunity for non-life-supporting devices

QEVO and QEVO ECU are intended for operation in an electromagnetic environment as specified below. The customer or the user of QEVO and QEVO ECU is responsible for ensuring that both products are used in such an environment.			
Electromagnetic immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment - guidelines
			Portable and mobile radio communication equipment should not be used closer to the QEVO and QEVO ECU including cables, than the recom-

			mended safety distance that is calculated using the equation applicable to the transmission frequency involved. Recommended safety distance
Conducted RF disturbances as per IEC 61000-4-6	3 V effective value 150 kHz to 80 MHz	3 V	$d = 1.2 \sqrt{P}$
Emission of RF disturbances in accordance with IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ for 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ for 800 MHz to 2.5 GHz where P is the output power rating of the transmitter in watts (W) according to the transmitter manufacturer's specifications and d is the recommended safety distance in meters (m). The field strength of stationary transmitters as determined by a site survey^a should be less than the compliance level in all frequency ranges.^b Interference may occur in the vicinity of equipment marked with the following symbol. 
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is influenced by absorption and reflection by structures, objects and persons.			
^a Theoretically, field strengths of stationary transmitters such as base stations for mobile telephones and mobile land radio equipment, amateur radio stations, AM and FM radio broadcast and TV broadcast transmitters cannot be predicted accurately. To assess the electromagnetic environment with respect to stationary RF transmitters, a site study of the electromagnetic phenomena should be considered. If the measured field strength in the location where the QEVO and QEVO ECU are used exceeds the compliance levels indicated above, the QEVO and QEVO ECU should be monitored to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the QEVO and QEVO ECU. ^b Field strengths should be less than 3 V/m over the frequency range from 150 kHz to 80 MHz.			

9.7.1.4 Recommended safety distances between portable and mobile RF communication equipment, QEVO and QEVO ECU

QEVO and QEVO ECU are intended for use in an electromagnetic environment in which RF disturbances are controlled. The customer or the user of the QEVO and QEVO ECU can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the QEVO and QEVO ECU - depending on the output power of the communication equipment as specified below.

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9.7.2 EMC Electromagnetic compatibility IEC 60601-1-2: 2014 (4th Edition)

The product and the control unit are subject to specific precautions with regard to electromagnetic compatibility (EMC) in the area of Professional Healthcare Facility Environment. In order to avoid the occurrence of EMC interference, the product may only be installed, operated and maintained in the manner indicated in these Instructions for Use and only with components supplied by ZEISS.

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- ▶ Only use accessories, transformers, cables and spare parts which are specified in these Instructions for Use or which are approved by ZEISS for this device.
- ▶ The product can be disturbed by other devices, also when these other devices conform with the emission requirements valid for them according to CISPR.
Do not operate the product in the direct vicinity of another device. If this cannot be avoided, you must monitor the proper function of the product and the other device.
- ▶ Do not use portable and mobile radios in the vicinity of the product. It cannot be ruled out that the function of the product will be impaired or that the performance of the devices deteriorates.
- ▶ Please note the EMC guidelines on the following pages.

9.7.2.1 Electromagnetic interference

QEVO and QEVO ECU are intended for operation in an electromagnetic environment as specified below. The customer or the user of QEVO and QEVO ECU is responsible for ensuring that both products are operated in such an environment.	
Interference measurements	Compliance
RF emissions as per CISPR 11	Group 1
RF emissions as per CISPR 11	Class A
Emission of voltage fluctuations/flicker as per IEC 61000-3-3	Not applicable
Harmonic emissions as per IEC 61000-3-2	Not applicable

NOTE

The emission properties of QEVO and QEVO ECU are suitable for use in industrial areas and hospitals (CISPR 11 Class A). If QEVO and QEVO ECU are used in a residential area (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting QEVO and QEVO ECU.

9.7.2.2 Electromagnetic immunity for all ME equipment and ME systems

QEVO and QEVO ECU are intended for operation in an electromagnetic environment as specified below. The customer or the user of QEVO and QEVO ECU is responsible for ensuring that both products are operated in such an environment.		
Electromagnetic immunity tests	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) as per IEC 61000-4-2	±8 kV contact discharge ±15 kV air discharge	±8 kV contact discharge ±15 kV air discharge
Fast transient/burst immunity as per IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines
Surges as per IEC 61000-4-5	±1 kV voltage Phase neutral conductor ±2 kV voltage Phase neutral conductor-earth	±1 kV voltage Phase neutral conductor ±2 kV voltage Phase neutral conductor-earth

Magnetic field for supply frequency (50/60 Hz) as per IEC 61000-4-8	30 A/m	30 A/m
Voltage dips, short interruptions and voltage variations as per IEC 61000-4-11	0 % U _T for 1/2 cycle	0 % U _T for 1/2 cycle
	0 % U _T for 1 cycle	0 % U _T for 1 cycle
	70 % U _T for 25/30 cycles	70 % U _T for 25/30 cycles
	0 % U _T for 250/300 cycles	0 % U _T for 250/300 cycles

9.7.2.3 Electromagnetic immunity for non-life-supporting devices

QEVO and QEVO ECU are intended for operation in an electromagnetic environment as specified below. The customer or the user of QEVO and QEVO ECU is responsible for ensuring that both products are used in such an environment.		
Electromagnetic immunity tests	IEC 60601 test level	Compliance level
Conducted RF disturbances as per IEC 61000-4-6	3 V 150 kHz to 80 MHz	3 V
	6 V ISM bands between 150 kHz and 80 MHz	6 V
Emission of RF disturbances in accordance with IEC 61000-4-6	3 V/m 80 MHz – 2.7 GHz	3 V/m
Radiated RF disturbance variables due to magnetic proximity fields of wireless communication equipment in accordance with EN 61000-4-3	27 V/m 385 MHz	27 V/m
	28 V/m 450 MHz, 810 MHz – 2.45 GHz	28 V/m
	9 V/m 710 MHz – 780 MHz, 5.24 GHz – 5.785 GHz	9 V/m

Empty page, for your notes

10 Accessories and spare parts

Ordering Data

Use only accessories and spare parts approved by ZEISS for this product. The safety of the product during operation cannot be guaranteed if accessories and spare parts which have not been approved by ZEISS are used.

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* Not available worldwide. For more information, please contact your local ZEISS sales representative.

Empty page, for your notes

11 Disposal

11.1 Disposal of the original packing

- Keep the original packing in the event of a relocation, maintenance, or repair.
- If you want to dispose of the original packing: Dispose of packing material by sending it for recycling through an acknowledged collection system.

11.2 Disposal of the product and control unit

In accordance with applicable EU guidelines and national regulations at the time at which the product was brought onto the market, the product specified on the consignment note is not to be disposed of via the domestic waste disposal system or municipal waste disposal facilities.

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If products which have not been reprocessed are returned to ZEISS Service, pathogenic germs may spread.
Cross-contamination can occur.

- Reprocess the product completely before sending it to ZEISS Service for disposal.

NOTE**Environmental damage due to improper disposal!**

If the product and control unit are improperly disposed of via domestic waste or via municipal waste disposal services, pollutants will enter the environment.

- ▶ If you wish to dispose of the product, return it to your local ZEISS Service.
 - ▶ If you wish to dispose of the control unit, inform your local ZEISS Service.
-
- ▶ Get in touch with your local ZEISS Service and inform them that you want to dispose of your product.
The procedure is based on the same principle as sending in the product for "Direct maintenance" [43].
 - ▶ If you resell the product or its components, you must inform the buyer that the product needs to be disposed of in accordance with the currently applicable regulations.
 - ▶ Get in touch with your local ZEISS Service and inform them that you wish to dispose of the control unit. The control unit is permanently installed in the KINEVO 900, and may only be removed and disposed of by ZEISS Service.

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IEC

International Electrotechnical Commission

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Specifications subject to change

G-30-1976-en - 8.0 - 2020-05-05

The Instruction for use «Special endoscopic camera QEVO and Special endoscopic camera control unit QEVO ECU» in accordance with national and international requirements is standard for all countries. In accordance with national requirements and in order to comply with Russian legislation the Instruction for use is supplemented by these Addendums, which are the parts of the instruction for use for the products supplied to the Russian Federation.

Addendum I to the document Instruction for use « Special endoscopic camera QEVO and Special endoscopic camera control unit QEVO ECU»

Special endoscopic camera QEVO and Special endoscopic camera control unit QEVO ECU
Classification According to IEC 60601

Device and control unit classification according to IEC 60601-1. QEVO and QEVO ECU belong to the following classes:

- The QEVO ECU has a class 1 protection against electric shock
- QEVO meets the requirements for type BF for protection against electric shock
- QEVO mode: short
- QEVO ECU operation mode: continuous.
- ingress protection (against water or solid particles): not applicable
- inappropriate for use in environments with high oxygen content

Device and control unit classification according to IEC 60601-2. QEVO and QEVO ECU belong to the following classes:

- Electromagnetic compatibility (EMC): complies with IEC 60601-1-2, Class A (according to CISPR 11)

Paragraph 9.1. "Regulatory Information" sub-paragraph "Classification of the Product and the Control Unit according to IEC 60601-1" should be supplemented with the following information:

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- ingress protection (against water or solid particles): not applicable
- inappropriate for use in environments with high oxygen content

In Instruction for use «Special endoscopic camera QEVO and Special endoscopic camera control unit QEVO ECU» mode of operation instead of "Intermittent" should read as

"Short"

In Instruction for use «Special endoscopic camera QEVO and Special endoscopic camera control unit QEVO ECU» mode of operation should read as

"Continuous"

Paragraph 9.1. «Ambient requirements for operation» should be supplemented with the following information:

Products that are operated inside human organs and cavities must be nondefective during operation at nominal temperature from 32 to 42°C.

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Disposal in accordance with the sanitary rules and regulations of the Russian Federation. Medical waste hazard class: Б.

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
Printed on 07. 06. 2011

Key to symbols

Different symbols used in this user's manual draw your attention to safety aspects and useful tips. The symbols are explained in the following.

Warning!

The **warning triangle** indicates potential sources of danger which may constitute a risk of injury for the user or a health hazard.

Caution:

The **square** indicates situations which may lead to malfunction, defects, collision or damage of the instrument.

Note:

The **hand** indicates hints on the use of the instrument or other tips for the user.

OPMI®

OPMI® is a registered trademark of Carl Zeiss.

Contents

– Key to symbols	2
Safety	5
– Directives and standards	5
– Notes on installation and use	7
– Warning labels and notes	8
Description	9
– Intended use	9
– Controls and connections	10
Preparations for use	14
– Mounting the mouth switch to the 180° tiltable tube	14
– Attachment of the mouth switch to the straight tube	16
Checklist	18
Maintenance/Miscellaneous	19
– Care of the unit	19
– Ordering data	20
– Accessories	20
– Ambient requirements	21
Index	23

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Issue 2.2
Printed on 07. 06. 2011

Safety

The instrument described in this manual has been developed and tested in accordance with Carl Zeiss safety standards and with national and international regulations. A high degree of instrument safety is thus ensured.

We would like to inform you on the safety aspects involved in operating the instrument. This chapter contains a summary of the most important precautions to be observed.

Further safety notes are also contained in other parts of this user's manual; they are marked with a warning triangle containing an exclamation mark as shown here. Please pay special attention to these safety notes.

Safety is only ensured when this instrument is operated properly. Please read through this manual carefully before turning the instrument on. Also read through the user's manuals of the other equipment used with this instrument. You may obtain further information from our service organization or authorized representatives.

Directives and standards

The unit described in the user manual has been designed in compliance with:

- EN
- IEC
- UL
- CSA

In accordance with Directive 93/42/EEC, Annex II, Article 3, the complete quality management system of the company Carl Zeiss has been certified by the DQS Deutsche Gesellschaft zur Zertifizierung von Managementsystemen mbH, a notified body, under registration number 250881 MP21.

- The instrument must be connected to a special emergency backup line supply in accordance with the regulations or directives which apply in your country.
- This is a class I instrument as defined by Directive 93/42 /EEC.

- Please observe all applicable accident prevention regulations.

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Notes on installation and use

Safe working order

- Do not operate the equipment contained in the delivery package in
 - explosion-risk areas,
 - the presence of inflammable anesthetics or volatile solvents such as alcohol, benzine or similar chemicals.
- Do not station or use the instrument in damp rooms. Do not expose the instrument to water splashes, dripping water or sprayed water.
- Do not force cable connections. If the male and female parts do not readily connect, make sure that they are appropriate for one another. If any of the connectors are damaged, have our service representative repair them.
- Do not use a mobile phone in the vicinity of the equipment because the radio interference can cause the equipment to malfunction. The effects of radio interference on medical equipment depend on a number of various factors and are therefore entirely unforeseeable.
- Modifications and repairs on these instruments or instruments used with them may only be performed by our service representative or by other authorized persons.
- The manufacturer will not accept any liability for damage caused by unauthorized persons tampering with the instrument; this will also forfeit any rights to claim under warranty.
- Use this instrument only for the applications described.
- This instrument is a high-grade technological product. To ensure optimum performance and safe working order of the instrument, its safety must be checked once every 12 months. We recommend having this check performed by our service representative as part of regular maintenance work.

If a failure occurs which you cannot correct using the trouble-shooting table, attach a sign to the instrument stating it is out of order and contact our service representative.


Warning labels and notes

Caution:

Observe all warning labels and notes!

If any label is missing on your instrument or has become illegible, please contact us or one of our authorized representatives. We will supply the missing labels.

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Description

Intended use

Pressing together the two parts of the mouth switch will release the magnetic brakes of the suspension system and enable you to guide the balanced surgical microscope almost effortlessly.

Warning!

For safety reasons, the system with mouth switch may only be used when it is correctly balanced!

With an incorrectly balanced system, brake release may lead to uncontrolled movements of the stand.

We would like to emphasize that the admissible total weight of the surgical microscope including accessories must not be exceeded when the mouth switch is used.

The accessory is intended for use with Carl Zeiss microscopes* in hospitals, clinics or other human medicine institutions.

The accessory must only be operated by physicians, nurses and other OR staff who have undergone appropriate training and observe the instructions of the user's manual. The installation conditions and the use of the system must meet microsurgical requirements:

- low vibration
- clean environment
- avoidance of extreme mechanical stress.

For more information, get in touch with our authorized representatives.


* OPMI Neuro, OPMI Neuro Multivision and OPMI Vario

Controls and connections

Mouth switch for 180° tiltable tube

- 1 Mouth switch
Pressing together the two parts of the mouth switch will release the magnetic brakes of the suspension system.
- 2 Rocker switch with click-stop
Press rocker switch downward - the mouth switch can then be swung to the side.
- 3 Clamping device for B-axis
to set the lower mouth switch tilt as required.
This screw can be removed completely to attach the mouth switch (1) to the holder for the straight tube (see page 16).
- 4 Clamping device for vertical adjustment
to set the height of the mouth switch as required.
- 5 Connector
for connection to the surgical microscope
- 6 Bores for fastening screws (3x)
for attachment to the 180° tiltable tube
- 7 Clamping device for A-axis
to adjust the tilt of the mouthpiece


Mouth switch holder for straight tube

- 1** Holding ring for straight tube
to mount the mouth switch to a straight tube, e.g. between straight tube and surgical microscope.
- 2** Toothed shaft
to vertically adjust the mouth switch as required
- 3** Clamping device for A-axis
to set the mouth switch distance as required.
This screw can be removed completely to attach the mouth switch to the toothed shaft.

Preparations for use


Mounting the mouth switch to the 180° tiltable tube

Note:

If the mouth switch has been mounted, the handgrip with asepsis cannot be attached to the right side of the 180° tiltable tube.

- Remove the three protection caps (2).
- Fasten the mouth switch holder (1) using three screws (3) contained in the delivery package.
- Tighten the screws using a 4 mm Allen key.
- Connect contact plug (4) to the socket for mouth switch (5) on the surgical microscope.
- Set height, tilt and distance of the mouth switch via the relevant clamping screws and properly tighten these screws.
- Make sure that the surgical microscope and the suspension system are properly balanced.

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Issue 2.2
Printed on 07. 06. 2011

Attachment of the mouth switch to the straight tube


- Remove clamping screw (5) from the toothed shaft (2).
- Loosen locking screw (4) by a few turns.
- Remove the mouth switch from the toothed shaft of the tiltable tube holder (delivery status).
- Align the mouth switch (3) and push it onto the toothed shaft (2).
- Clamp the mouth switch (3) using clamping screw (4).
- Insert the clamping screw (5) into the toothed shaft (2) again and clamp it
- Mount the holding ring (6) to the front flange of the straight tube.
- Make sure that the two recesses (8 and 9) coincide.
- Place straight tube with complete mouth switch on the microscope body and firmly tighten securing screw (10).

Note:

You can install further accessories between the straight tube and the microscope body. Lock these units in position in the same way using securing screw (10).

- Connect contact plug (7) to the socket for mouth switch on the surgical microscope.

Checklist

Always check the following points before use (without patient!):

- The 180° tiltable tube or straight tube is securely seated.
- The mouth switch holder has been securely mounted to the respective tube.
- The contact plug of the mouth switch has been properly connected to the socket for mouth switch on the surgical microscope.
- Make sure that the mouth switch is in a position convenient for you.
- Make sure that there is enough room between the instrument and the accessories installed.
- Turn on the instrument at the power switch of the suspension system.
- Check that the surgical microscope has been correctly balanced.

Warning!

For safety reasons, the system must only be used when correctly balanced!

With an incorrectly balanced system, brake release may lead to uncontrolled movements of the stand.

For this reason, hold the surgical microscope securely at its handgrips before releasing the magnetic brakes.

- Check whether the mouth switch functions properly, i.e. all magnetic brakes of the system are released when the mouth switch is pressed together and all magnetic brakes of the system are locked when the mouth switch is released.

Warning!

If the mouth switch fails, you must not use the instrument for safety reasons. Repair the cause of the defect or call our service department.

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Ordering data

Mouth switch

Description	Cat. No.
Mouth switch for 180° tiltable tube	1177-805
Mouth switch holder for straight tube (Requires mouth switch 1177- 805 to be ordered!)	1116-378

Accessories

Please observe the following:

Only operate the instrument with the accessories included in the delivery package. If you want to use other accessories, make sure that Carl Zeiss or the manufacturer of the accessories has proved and confirmed that these accessories meet the respective technical safety standards and can be used without risk.

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Issue 2.2
Printed on 07. 06. 2011

Ambient requirements

For operation	Temperature Rel. humidity Air pressure	+10 °C...+40 °C 30%...75% 700 hPa...1,060 hPa
For transportation and storage	Temperature Rel. humidity (without condensation) Air pressure	- 40 °C...+70 °C 10%...100% 500 hPa...1,060 hPa

The unit meets the essential requirements stipulated in Annex I to the 93/42/EEC Directive governing medical devices. The unit is marked with:

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Index

A

Accessories 20
Accident prevention regulations 6
Ambient requirements 21
Attachment of mouth switch to the straight tube 16

C

Care of the unit 19
CE marking 21
Checklist 18
Controls and connections 10

D

Directive 93/42/EEC 21
Directives 5

E

Emergency backup line supply 5

I

Installation, use 7
Instrument class 5

K

Key to symbols 2

M

Mechanical surfaces, cleaning 19
Mobile phones 7
Mounting the mouth switch to the 180° tiltable tube 14

N

Notes on installation 7
Notes on safety 5
Notes on use 7

O

Ordering data 20

P

Preparations for use 14

Q

Quality management system 5

S

Safety 5
Safety check 7
Safety standards 5

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FOR THE RUSSIAN FEDERATION ONLY

The Instruction for use «Mouthswitch for 180° Tilttable tube» in accordance with national and international requirements is standard for all countries. In accordance with national requirements and in order to comply with Russian legislation the Instruction for use is supplemented by these Addendums, which are the parts of the instruction for use for the products supplied to the Russian Federation.

For Russian Federation only one variant is available: Mouthswitch for 180° Tilttable tube

Addendum I to the document Instruction for use «Mouthswitch for 180° Tilttable tube»

Paragraph «Care of the unit» (page 19) should be supplemented with the following information:

All mechanical surfaces of the system can be cleaned by wiping them with a damp cloth. Do not use any aggressive or abrasive cleaning agents. Remove any possible residue using a mixture of 50% ethyl alcohol and 50% distilled water plus a dash of household dishwashing liquid.

Clearing: Before the first use and after each use.

Paragraph «Ambient requirements» in disposal part should be supplemented with the following information:

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Instruction for Use (IFU)

Dear Madam or Sir,

Please find enclosed the English and Russian Instruction for Use (IFU) for INFRARED 800 with FLOW 800 Option.

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Address of Record:
Goeschitzer Strasse 51 - 52
07745 Jena, Germany

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Aalen, den 04.05.2021

Kurz, Notar

INFRARED 800 with FLOW 800 Option For KINEVO 900

Instructions for Use

Copyright

© 2020, Carl Zeiss Meditec AG, Jena

Brand names

KINEVO 900 and the fluorescence application INFRARED 800 with FLOW 800 Option are trademarks or registered trademarks of Carl Zeiss Meditec AG or another company of the ZEISS group in Germany and/or other countries.

All other trademarks mentioned in this document are the property of their respective owners.

Table of contents

1	Notes on the instructions for use	5
1.1	Product name	5
1.2	Area of application	5
1.3	Purpose and storage of the documentation	6
1.4	Questions and comments	6
1.5	Conventions in this document	6
1.5.1	Conventions in all text areas	6
1.5.2	Conventions in a course of action	7
1.6	Other applicable documents	7
2	Safety notes	9
2.1	Target group	9
2.2	Area of use	9
2.2.1	Intended use	9
2.2.2	Normal use	9
3	Product description	11
3.1	Labeling on KINEVO 900	11
3.2	Functional description	12
3.3	Licenses	15
3.4	Software description	16
3.4.1	KINEVO 900 user interface	16
3.4.2	"INFRARED 800" menu	18
3.4.3	"FLOW 800" menu	20
4	Configuring INFRARED 800	27
4.1	Configuring the [Fluorescence] button on the handgrip and (optionally) on the foot control panel	27
4.2	Adjusting the screen display	28
4.3	Setting the fluorescence parameters	28
4.4	Showing INFRARED 800 with FLOW 800 (option) in the eyepieces	31
5	Prior to each use	33
5.1	Performing a function test	33
5.2	Preparing patients	34
6	Operation	35
6.1	Operating INFRARED 800 with FLOW 800 Option via handgrip or foot control panel	37
6.2	Operating INFRARED 800 with FLOW 800 Option via the monitor	39

6.3	Analyzing and editing FLOW 800 data	42
6.3.1	Calling FLOW 800 data.....	42
6.3.2	Observing different characteristics of the FLOW 800 display.....	43
6.3.3	Setting markers (ROIs)	45
6.3.4	Interpreting and editing the diagram.....	47
6.3.5	Optimize View.....	50
6.3.6	Comparing FLOW 800 displays.....	52
6.3.7	Transferring markers (ROIs).....	54
7	Technical Data KINEVO 900	55
7.1	Integrated HD camera	55
7.2	Digital video outputs	55
7.3	4K video outputs	55
7.4	Digital video inputs.....	55
7.5	Ambient requirements for operation	56
7.6	Ambient requirements for transport and storage.....	56
7.7	Operation with intraoperative MRT systems	56
8	Order data.....	57
8.1	Fluorescence target	57
8.2	Licenses	57
	Glossary	59
	Keyword index	61

1 Notes on the instructions for use

1.1 Product name

INFRARED 800 is an application for KINEVO 900.

FLOW 800 is an optional upgrade for the application INFRARED 800.

All applications and upgrades are activated by their respective licenses on the KINEVO 900.

Display and operation

Only activated applications can be displayed and operated on the KINEVO 900 monitor.

1.2 Area of application

These instructions for use apply to INFRARED 800 with FLOW 800 Option and the following identification label:

- reference number: 7012

These instructions for use describe only the INFRARED 800 with FLOW 800 Option fluorescence application. Read and follow the applicable instructions for use for KINEVO 900.

The INFRARED 800 with FLOW 800 Option fluorescence technologies are fully integrated in the KINEVO 900 and are available for sale in the United States. They are also available in other markets.

Depending on the national regulations of the respective country where the system is operated, an additional approval may be required for the system and its application. Please contact your local ZEISS contact person to obtain more information regarding the registration status of INFRARED 800 with FLOW 800 Option in your country.

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1.4 Questions and comments

- ▶ If you have any questions or comments concerning these instructions for use or the device itself, please contact ZEISS Service.

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This is a cross-reference: Questions and comments [6].

This is **bold type**.

This is a `software code or program text`.

Names of software dialogs, fields or menus, and software messages, are marked by quotation marks:

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☒ This is a requirement that must be met before the start of a sequence of actions.

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1.6 Other applicable documents

Please observe the applicable instructions for use for KINEVO 900 in the software version used.

Empty page, for your notes

2 Safety notes

2.1 Target group

The present Instructions for Use are intended for physicians, medical and technical staff and nurses who prepare and operate the product after receiving appropriate instruction.

It is the duty of the equipment owner/operator to train and instruct all operating personnel.

The intended usage site for the product is operating rooms in clinics and hospitals.

2.2 Area of use

2.2.1 Intended use

INFRARED 800 with FLOW 800 Option used during fluorescence guided surgery aids in viewing the visual assessment of intra-operative blood flow as well as vessel patency in surgical procedures in neurosurgery, plastic and reconstructive procedures and coronary artery bypass graft surgery.

2.2.2 Normal use

INFRARED 800 with FLOW 800 Option is an accessory product for a surgical microscope that is used to observe and visually assess the intraoperative blood flow in the region of the cerebral blood vessels. Its applications include i.a. the assessment of cerebral aneurysms and vascular occlusions as well as the permeability of very small perforator vessels. In addition, it supports real-time visualization of the blood flow and visual assessment of vessel types before and after surgical procedures for arteriovenous malformations (AVM).

The use of INFRARED 800 with FLOW 800 Option during fluorescence-supported surgery also assists in the viewing and assessment of the interoperative blood flow and of the vessel patency for surgical bypass procedures in neurosurgery, plastic and reconstructive surgery, and coronary artery bypass graft surgery.

When using INFRARED 800 in conjunction with a fluorescent agent, the medical contraindications which must be considered when using suitable fluorescent agents and state-of-the-art examination techniques apply.

- Check whether the fluorescent agent can be excited in the wavelength range from 700 to 780 nm and whether it emits fluorescent light of sufficient intensity in the wavelength range from 820 to 900 nm.

The automatic detection (AD) and the FLOW 800 display require a correspondingly high fluorescence signal. If the fluorescence signal is too low, this can lead to an inaccurate AD calculation and FLOW 800 display.

As in almost all procedures, false-positive and false-negative results may occur. Evaluation using other methods is necessary.

3 Product description

3.1 Labeling on KINEVO 900

The labeling for INFRARED 800 with FLOW 800 Option is attached to the KINEVO 900.


Figure 1: Labeling for INFRARED 800 with FLOW 800 Option on the KINEVO 900

Pos.	Symbol	Designation	Explanation
1		Identification label INFRARED 800	Provides information on: ■ Product designation ■ Manufacturer ■ Item and serial number ■ Date of manufacture ■ UDI labeling ■ CE label
		Serial No.	-
		Catalog number	Reference number:
		Date of manufacture	Contains the product's date of manufacture.
		UDI	UDI labeling

3.2 Functional description

The INFRARED 800 fluorescence application supports the observation and visual assessment of the intraoperative blood flow as well as the vascular permeability during surgical procedures. This process records and saves video data which can be played back during or after the surgical procedure and forms the basis for a detailed analysis of the fluorescence application INFRARED 800 with the FLOW 800 module.

The FLOW 800 module summarizes the video data recorded by the INFRARED 800 fluorescence application. This summary of the video data supports fast and objective interpretation of the blood flow. Markers (ROIs) can be set at important and/or interesting locations in the tissue. This allows the different characteristics of the influx and flow properties of the blood at these marked locations to be conveniently interpreted.

Another important function is the Compare Mode of the FLOW 800 module. In the Compare Mode, the currently recorded video data can be compared with a patient's previously recorded video data. The Compare Mode can be used, for example, during an AVM procedure to visually assess vessel types and their blood flow and compare them at the beginning of and after surgery.

Composition of the application

The application INFRARED 800 with FLOW 800 Option comprises two modules:

- INFRARED 800 provides the hardware and software for the video recording of the infrared fluorescence during the operation.
- The optional module FLOW 800 provides an easy-to-use interface that enables comfortable visualization of the fluorescence video data recorded by the INFRARED 800 module. It compiles the INFRARED 800 video data in color coded overviews in which the sequence in which the fluorescence appears ("Delay" characteristic), the intensity of the fluorescence ("Intensity" characteristic) and the speed of the fluorescence rise ("Speed" characteristic) are displayed in a convenient and comprehensible manner.

Both modules are integrated in the KINEVO 900. You can open and operate both modules via the handgrips, the foot control panel or the main monitor of the KINEVO 900. First, you record video data during the infrared fluorescence with the INFRARED 800 module. If necessary, you can use the "Replay" function for further assessment of the blood flow and, optionally, FLOW 800 for further interpretation.

Main functions

Designation	Functional description
INFRARED 800	Records video data of the fluorescing tissue in the surgical field during the operation. The video recording process is divided into three phases: ■ Setup - Preparation phase; here you can set the focus and zoom parameters (provided that the "Skip Setup" function has been deactivated) and select the image section; then the recording phase starts. ■ Record - Recording phase; starts the detection and display of the fluorescence signals and saves this as a video file. ■ Replay - Playback phase; plays the recorded video files on the screen: "Short replay" plays a short sequence at the beginning of the video (duration adjustable) on repeat (repetition rate adjustable). "Long replay" plays the entire video on repeat (repetition rate adjustable). If the KINEVO 900 is equipped with the 3D option, you can simultaneously record a white light video with the INFRARED 800 module and record a INFRARED 800 video and display both recordings on the screen as an overlay or picture-in-picture.
Automatic Gain Control	Automatic camera gain increase during the recording phase. You can switch this function on or off.
AutoDetection	The system determines the point of time in the video data when the fluorescence agent flows into the tissue. If this function is switched on, the playback phase automatically starts at this point of time. Video data which precede this point of time are not displayed. You can switch this function on or off.
AutoZoom	The system starts the recording phase at the magnification previously set in the setup menu. You can switch this function on or off.
Analytical characteristics of the FLOW 800 module	The FLOW 800 module processes the video data of the fluorescing tissue according to the following characteristics: ■ Delay of the fluorescence influx into the tissue ■ Maximum intensity achieved by the fluorescence signal in the tissue ■ Speed at which the fluorescence signal increases (supports the assessment of the speed of blood flow)
ROI (region of interest)	During the playback of the video data, you can set markers (ROIs: regions of interest) at locations which are of interest to you. The intensity of the fluorescence in the tissue is processed and displayed as a curve in a diagram for each marker set. With the help of this diagram, you can conveniently and quickly evaluate and compare the time course of the fluorescence intensity in the tissue.

Designation	Functional description
Diagram display	The diagram displays the time course of the fluorescence intensity for the previously set markers (ROIs). You can observe the change of intensity of the fluorescence brightness (in %) over time (in sec.) in the diagram. Each marker (ROI) is displayed as a separate curve in the diagram. In addition, the following values are displayed: ■ "Delay" characteristic ■ "Speed of fluorescence increase" characteristic ■ Time to Peak (TTP) (Period of time the contrast medium requires to reach its maximum concentration in the tissue (time until maximum hyperintensity)) ■ Rise Time (RT) (Period of time the contrast medium requires until an increase from a preset start value to a preset maximum value occurs)
Auto-correction	This function compensates for small movements of the KINEVO 900 (e.g. vibrations) and, if necessary, for poor illumination towards the edges of the image (vignetting).
Optimize View	Using this function, you can select a partial sequence of the INFRARED 800 video recording that is to be summarized in the FLOW 800 module. In this way you can e.g. exclude artifacts which occurred at the beginning or at the end of the INFRARED 800 video recording from the summary, thus optimizing the display. You then can save the optimized display; the entire video recording also will be preserved.
Compare FLOW 800 displays	Using this function, you can compare two FLOW 800 displays of the same patient recorded at different times. Thus, for example, you can record and save FLOW 800 displays of your patient before and after surgery to enable a before/after comparison of the fluorescence intensity in the tissue. In the compare view, you can display the fluorescence intensity based on the Delay, Intensity and Speed characteristics. You can compare different diagrams only if you have saved a display of each of the diagrams to be compared.
Transfer markers (ROIs)	Using this function, you can transfer the markers of a previously recorded FLOW 800 display to your currently recorded FLOW 800 display. This enables you to conveniently compare the fluorescence intensity of the previous FLOW 800 display with that of the current one based on the Delay, Intensity and Speed characteristics. For example, you can compare the fluorescence influx into the tissue before and after surgery side-by-side.

3.3 Licenses

Each module requires a license which is managed by KINEVO 900. You can open the modules and operate them only if you have valid licenses for this purpose:

Designation	Symbol	Scope of operation
INFRARED 800 Compact		Recording of video data during the operation without the main functions AutoGain and AutoDetection, i.e. the camera gain and recording time point must be set manually. The "Replay" playback phase is not available.
INFRARED 800		Recording and saving of video data during the operation with the main functions AutoGain and AutoDetection [12].
INFRARED 800 with FLOW 800		Playback, analysis and editing of the recorded video data.

To determine which licenses are activated on your KINEVO 900, check the identification label for applications [11] on the KINEVO 900 or open the "Fluorescence" function on the user interface .

3.4 Software description

3.4.1 KINEVO 900 user interface

The INFRARED 800 with FLOW 800 Option application can be opened, controlled and configured by using the following control elements on the KINEVO 900:

- Main monitor (touchscreen): Open, control and configure INFRARED 800 with FLOW 800 Option application
- Handgrips: Open and control INFRARED 800 with FLOW 800 Option application
- Foot control panel: Open and control INFRARED 800 with FLOW 800 Option application

The user interface is subdivided into different areas with different controls.


Figure 2: Overview of user interface and status bar

1	Live video camera image	2	Status bar
3	Home (Navigation for video camera live image)	4	Edit patient data
5	Display of patient name, opens patient archive	6	Edit user data
7	Display of user name, opens user archive	8	Status information (opens corresponding menu)
9	Date, time	10	Setup menu
11	Help	12	Extended operation (open)

13	Main menu bar	14	Start INFRARED 800 with FLOW 800 Option application
15	Number of images taken	16	Capture still image (photo)
17	Number of videos recorded	18	Video recording Start/Stop

All control elements, the status bar and the main menu bar can be hidden by tapping on the live image (full screen mode). Only the video runtime is displayed during a recording in full screen mode. If you tap on the monitor again, all control elements, the status bar and the main menu bar will reappear.

The displayed values are rounded. They serve only to enhance the display and do not represent a measurement function.

Detailed information on the user interface is provided in the instructions for use of KINEVO 900 in the software version used.

3.4.2 "INFRARED 800" menu

When the INFRARED 800 application is started, the "INFRARED 800" menu opens. The top edge of the screen shows the different phases of the INFRARED 800 application. The blue bar marks the phase currently in progress.


Figure 3: "INFRARED 800" menu

Pos.	Name	Explanation
1	"Setup" preparation phase	In the "Setup" preparation phase you can set the focus and zoom parameters and select the image section.
2	"IR Record" recording phase	During the recording phase, the fluorescence signal is detected, displayed and recorded as a video file. During the recording phase, the current values for the light intensity and camera gain are displayed.
3	"Short Replay" playback phase	The "Short Replay" playback mode repeatedly and successively plays back a short sequence at the beginning of the video (adjustable duration and adjustable repetition rates).
4	"Long Replay" playback phase	The "Long Replay" playback mode repeatedly and successively plays back the complete video (adjustable repetition rates).
5	Marking of the recording area	The two vertical white broken lines mark the area of the INFRARED 800 video in which summaries are created for the FLOW 800.

Pos.	Name	Explanation
6	Take photo of current display	Creates and saves documentation for the current display. You can view these images in the "Patient" menu.
7	Video recording Start/Stop	You can optionally start or stop the INFRARED 800 video recording by tapping this button.
8	Notification window	This window displays notifications concerning recommended settings that may help improve the INFRARED 800 recordings.
9	Zoom	Information indicator for the "Zoom" function. If the indicator line is colored white, the zoom is currently within the recommended zoom range.
10	Focus	Information indicator for the "Focus" function. If the indicator line is colored white, the focus is currently within the recommended focus range.

3.4.3 "FLOW 800" menu

Once a patient's video data have been recorded with INFRARED 800, you can display them in the FLOW 800 menu.


Figure 4: Structure of the "FLOW 800" menu

Pos.	Name	Explanation
1	Scale with slide control	Slider for individual adjustment of the color coding. The slider can be adjusted continuously.
2	Delay characteristic	Displays the time course (in sec.) of the fluorescence influx into the tissue. The red colored areas are perfused first; the blue colored areas are perfused later (delayed).
3	Intensity characteristic	Displays the maximum fluorescence achieved in the tissue. Values are shown as % of maximum value (100%). The fluorescence intensity is high in the red colored areas and low in the blue colored areas.
4	Speed characteristic	Displays the speed (AI/sec.) of the brightness rise in the tissue. The rise rate is high in red colored areas and low in blue colored areas.
5	Diagram	The diagram displays the time course of the fluorescence intensity for the previously set markers (ROIs). You can observe the change of intensity of the fluorescence brightness (in %) over time (in sec.) in the diagram. Each marker (ROI) is displayed as a separate curve in the diagram.

Pos.	Name	Explanation
6	Compare Mode	You can start the "Compare Mode" for any characteristic (Delay, Intensity, Speed). You can compare the diagram view of the current recording with the diagram views of previously created recordings of the same patient, provided that you previously saved the diagram views of the other recordings.
7	ROI (marker)	Tap on one or more points on the screen which are of interest to you and you would like to analyze. You can set a maximum of eight markers (ROIs). The following options are available in the "ROI" field: ■ Show/hide ROI ■ Delete ROI
8	Options for ROIs	After displaying the options in the "ROI" field, you can set the following parameters for each ROI: ■ Size ■ Shape (square or adaptive) ■ Color
9	Optimize View	Using this function, you can select a partial sequence of the INFRARED 800 video recording that is to be summarized in the FLOW 800 module. In this way you can e.g. exclude artifacts which occurred at the beginning or at the end of the INFRARED 800 video recording from the summary, thus optimizing the display. You then can save the optimized display; the entire video recording also will be preserved.
10	Auto-correction	This function compensates for and corrects small movements of the KINEVO 900 (e.g. vibrations) and, if necessary, for poor illumination towards the edges of the image (vignetting).
11	Take photo of current display	Creates and saves documentation for the current display. You can view these images in the "Patient" menu.

3.4.3.1 "Diagram" menu

The diagram displays the time course of the fluorescence intensity in the tissue for the previously set markers (ROIs). You can observe the change of fluorescence intensity (in %) over time (in sec.) in the diagram. Each previously set marker (ROI) is displayed as a separate curve.


Figure 5: Structure of the "Diagram" menu

Item	Name	Explanation
1	"Intensity" slide control	Slide control for scaling the fluorescence brightness (in %). You can set the limits of the display area between 0% and 100%.
2	Influx characteristics of the fluorescent agent	Select the curve of a marker (ROI) and read off the influx characteristics of the fluorescent agent in this region: ■ Delay ■ Speed ■ Time to Peak (period of time between AutoDetection and maximum achieved fluorescence intensity (time until maximum hyperintensity)) ■ Rise Time (period of time in which the fluorescence intensity increases from 10% to 90% of the maximum value achieved)
3	Additional "ROI" menu	The following options are available in the "ROI" field: ■ Show/hide ROI ■ Delete ROI

Item	Name	Explanation
4	Optimize View	Using this function, you can select a partial sequence of the INFRARED 800 video recording that is to be summarized in the FLOW 800 module. In this way you can e.g. exclude artifacts which occurred at the beginning or at the end of the INFRARED 800 video recording from the summary, thus optimizing the display. You then can save the optimized display; the entire video recording also will be preserved.
5	Auto-correction	This function compensates for small movements of the KINEVO 900 (e.g. vibrations) and, if necessary, for poor illumination towards the edges of the image (vignetting).
6	Video view	The INFRARED 800 video can be played back in the video view. The point of time currently displayed in the video is indicated by a white bar in the diagram. This enables you to compare the video display with the diagram display.
7	"Video playback" slide control	The slide control shows the current position in the video. You can adjust the slider in order to continuously fast-forward or rewind the video view.
8	Fast-forward/Rewind	You can move the video view either one sequence forward or one sequence back.
9	Start	To start the video view, tap on the [Start] button.
10	Take photo of current display	Creates and saves documentation for the current display. You can view these images in the "Patient" menu.
11	"Time" slide control	Time axis scaling. You can set the limits of the display area between the start and the end of the fluorescence recording.

3.4.3.2 “Optimize View” menu

Using this function, you can select a partial sequence of the INFRARED 800 video recording that is to be summarized in the FLOW 800 module. In this way you can e.g. exclude artifacts which occurred at the beginning or at the end of the INFRARED 800 video recording from the summary, thus optimizing the display. You then can save the optimized display; the entire, unedited video recording will nevertheless be preserved.


Figure 6: Structure of the “Optimize View” menu

Item	Name	Explanation
1	Video	The INFRARED 800 video displays the video recording.
2	Diagram	The diagram displays the time course of the fluorescence intensity for the previously set markers (ROIs). You can observe the change of intensity of the fluorescence brightness (in %) over time (in sec.) in the diagram. Each marker (ROI) is displayed as a separate curve in the diagram.
3	Fast-forward/Rewind	Navigate forward or in reverse frame by frame in the video sequence.
4	Cut out	Determines the end point of the timeline selection window, and therefore of the scheduled video recording.
5	Cut in	Determines the starting point of the timeline selection window, and therefore of the scheduled video recording/display.
6	Time marker	Marks the current position in the timeline selection window during video playback.

Item	Name	Explanation
7	"Video playback" slide control	The gray slider keeps moving throughout the entire length of the video recording. Slide the gray slider to jump to an interesting area in the video recording without using the [Fast-forward] / [Rewind] buttons.
8	AutoDetection (AD) mark	Marks the point of time when the contrast medium flows into the tissue.
9	Video playback time	Displays the current playback time of the video recording.
10	Duration of video sequence cutout	Total duration of the video sequence which you cut out.
11	Timeline selection window	Displays sequences of the video. The start and end of the timeline selection window is bounded by blue bars. You can slide these blue bars to select and cut the sequence relevant for you in the video. The FLOW 800 views are then recreated based on the image data contained in the sequence you selected.
12	Start	The video recording is played back several times in succession without any limitation of the playback time. Tapping the button again stops the playback of the video recording.

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User-specific settings

User-specific settings are marked by a green dot placed in front of their designation. These changes in settings are adopted only for the current user.

4.1 Configuring the [Fluorescence] button on the handgrip and (optionally) on the foot control panel

Prerequisite

- ☒ KINEVO 900 is switched on and ready for operation
- ☒ Optional: foot control panel is paired with KINEVO 900

Action

1. To assign the "Fluorescence" function to a handgrip button, tap on  Settings →  Operation.
2. Tap on button A or B (you can assign the "Fluorescence" function only to button A or button B on the handgrip).
 - ⇒ The available functions appear in the selection field on the right side.
3. Tap on the "Fluorescence" function in the selection field.
 - ⇒ The "Fluorescence" function is accepted for the desired button and displayed next to the button symbol.
4. To assign the "Fluorescence" function to a foot control panel button, tap on  Settings →  Operation.
5. Tap on a freely configurable button.
 - ⇒ The available functions appear in the selection field on the right side.
6. Tap on the "Fluorescence" function in the selection field.
 - ⇒ The "Fluorescence" function is accepted for the desired button and displayed next to the button symbol.

4.2 Adjusting the screen display

You can adjust the live display of your INFRARED 800 fluorescence application on the main screen. You can select the following display modes:

- Full Screen: Displays the fluorescence application on the entire main screen
 - Live Overlay (available only for 3D HD option): Displays the fluorescence application in the eyepieces of the KINEVO 900
 - Picture-in-Picture (available only for 3D HD option): White light video and INFRARED 800 video are displayed on the main screen at the same time
- ☒ KINEVO 900 is switched on and ready for operation

Prerequisite

Action

1. Tap on  Settings  Displays.
2. In the "INFRARED 800 Live Display" field, select the way you want to have your INFRARED 800 fluorescence application displayed on the main screen.
 - ⇒ You can choose between the "Full Screen", "Live Overlay" (available only for 3D HD option) and "Picture-in-Picture" (available only for 3D HD option) display modes.

4.3 Setting the fluorescence parameters

You can set the following parameters for the INFRARED 800 with FLOW 800 application:


- AutoGain: Automatic camera gain increase during the recording phase.
- AutoDetection: The system determines the point of time when the fluorescent agent flows into the tissue from the video data. If this function is switched on, the playback phase automatically starts at this point of time. Video data which precede this point of time are not displayed.
- Skip Setup: The preparation phase during which you can once again adapt the focus and zoom parameters is skipped.
- AutoZoom: The system starts the recording phase at the magnification previously set in the setup menu.
- Short Replay (number and duration): Repeatedly and successively plays back a short sequence at the beginning of the video (adjustable duration and adjustable repetition rates).
- Long Replay (number): Repeatedly and successively plays back the complete video (adjustable repetition rates).

Prerequisite

- ☒ INFRARED 800 with FLOW 800 license is available
- ☒ The KINEVO 900 is switched on and ready for operation


Action

1. Tap on  Settings  → Fluorescence.
2. Select INFRARED 800 or INFRARED 800 with FLOW 800 in the "Fluorescence" field (only the modules which are activated on your KINEVO 900 are displayed).
⇒ The button on the handgrip or foot control panel to which you assigned the "Fluorescence" function is now active.
3. To switch on the automatic camera gain, activate the "AutoGain" slide switch in the "IR 800" field.
⇒ The button is blue  : The camera gain is controlled automatically during the recording phase.
4. To switch off the automatic camera gain, deactivate the "AutoGain" slide switch in the "IR 800" field.
⇒ The button is gray  : The automatic camera gain is switched off and must be controlled manually.
5. To switch on the "AutoDetection" function, activate the "AutoDetection" slide switch in the "IR 800" field.

⇒ The button is blue  : AutoDetection is switched on. The black leader (X) without any fluorescence signal is hidden in the playback phase.

6. To switch off the "AutoDetection" function, deactivate the "AutoDetection" slide switch in the "IR 800" field.

⇒ The button is gray  : AutoDetection is switched off. The black leader (X) is also played back in the playback phase.

7. To skip the "Setup" phase: Activate the "Skip Setup" slide switch in the "IR 800" field.

⇒ The button is blue  : The INFRARED 800 video recording starts the first time the button is pressed.

8. To not skip the "Setup" phase: Deactivate the "Skip Setup" slide switch in the "IR 800" field.

⇒ The button is gray  : The first time you press the button, you activate a Setup phase in which you can set the image section, focus and zoom.

9. To switch on the "AutoZoom" function, activate the "AutoZoom" slide switch in the "IR 800" field.

⇒ The button is blue  : The total magnification is automatically changed to the preconfigured start value.

10. Adjust the start value of the automatic zoom with the "AutoZoom (start value)" control in the "IR 800" field.

11. To switch off the "AutoZoom" function, deactivate the "AutoZoom" slide switch in the "IR 800" field.

⇒ The button is gray  : A setup dialog for changing the total magnification outside of the preconfigured range is displayed.

12. Set the start value of the automatic and manual camera gain with the "Gain (start value)" control in the "IR 800" field.

13. Set the desired repetition rate for the short replays with the "Short Replay [Number]" control in the "IR 800" field.

⇒ Short replays are played back successively until the repetition rate set by you is reached or the Fluorescence button on the handgrip or foot control panel is pressed. You can set the repetition rate to a value ranging from 0 to 10.

14. Set the duration of the short replays with the "Short Replay [Duration]" control in the "IR 800" field.

15. Set the desired repetition rate for the long replays with the "Long Replay [Number]" control in the "IR 800" field.

⇒ Long replays are played back successively until the repetition rate set by you is reached or the Fluorescence button on the handgrip or foot control panel is pressed. You can set the repetition rate to a value ranging from 0 to 10.

4.4 Showing INFRARED 800 with FLOW 800 (option) in the eyepieces

You can show the INFRARED 800 and FLOW 800 modules in the eyepieces of the KINEVO 900 by using the "MultiVision" function.

Note:

If the recorded FLOW 800 displays have not yet been calculated at the end of the "Long Replay" playback mode, the shutter opens. In this case, the recorded FLOW 800 displays are not automatically injected into and shown in the eyepieces in order to ensure that the surgeon's vision is not impaired by a sudden and unexpected display.

Prerequisite

- ☒ A license for the "MultiVision" function is present
- ☒ The KINEVO 900 is switched on and ready for operation

Action

1. To assign the "MultiVision" function to a handgrip button, tap on  Settings →  Operation.
2. Tap on button A or B (you can assign the "MultiVision" function only to button A or button B on the handgrip).
 - ⇒ The available functions appear in the selection field on the right side.
3. Tap on the "MultiVision" function in the selection field.
 - ⇒ The "MultiVision" function is accepted for the desired button and displayed next to the button symbol.
4. To assign the "MultiVision" function to a foot control panel button, tap on  Settings →  Operation.
5. Tap on a freely configurable button.
 - ⇒ The available functions appear in the selection field on the right side.
6. Tap on the "MultiVision" function in the selection field.
 - ⇒ The "MultiVision" function is accepted for the desired button and displayed next to the button symbol.
7. Tap on  Settings →  MultiVision.
8. Switch on the "MultiVision" function by setting the slide switch in the "MultiVision" field to "On".
 - ⇒ The slide switch is blue  : The button on the handgrip or foot control panel to which you assigned the "Multi-Vision" function is now active.
9. Switch on the Live Overlay function for the INFRARED 800 fluorescence application by activating the slide switch in the "IR 800" field.

⇒ The slide switch is blue  : During the INFRARED 800 fluorescence application, you can see the live image of the surgical field in the eyepiece and the video recording in the INFRARED 800 mode overlaid above it.

10. Switch on the overlay of the FLOW 800 module by activating the slide switch in the "FLOW 800" field.

⇒ The slide switch is blue  : After the INFRARED 800 fluorescence application, you can see the FLOW 800 display in the eyepieces.

5 Prior to each use

NOTE

Malfunction due to failure of the light source or the electronics!

If the light source or the electronics in the KINEVO 900 fail, the INFRARED 800 with FLOW 800 Option fluorescence application cannot be completed.

- ▶ Prior to every operation and/or use of the INFRARED 800 with FLOW 800 Option fluorescence application, take suitable precautions to ensure that you will be able to complete the operation even if the INFRARED 800 with FLOW 800 Option fluorescence application fails.

NOTE

Malfunction due to additional infrared light sources!

Infrared light sources (e.g. navigation systems) may cause interfering signals during the INFRARED 800 fluorescence application! ZEISS recommends taking the following precautions to prevent this:

- ▶ Ensure that no infrared cameras of navigation systems are aimed at the surgical field while INFRARED 800 is being used.
- ▶ Before using INFRARED 800, check whether any light sources located in the room are causing interfering signals. If possible, position any such light sources further away from the surgical field.
- ▶ Prevent any reflections of infrared light from entering the eyepieces.

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The INFRARED 800 fluorescence target and the touchscreen are unsterile!

- ▶ Perform the function test prior to the operation (without a patient and outside of the sterile area!).

NOTE

Malfunction due to weak fluorescence signal!

The automatic display of the fluorescence rise (AutoDetection (AD)) and the FLOW 800 display require a clear fluorescence signal. If the fluorescence signal is too weak, this may lead to an inaccurate display.

- ▶ Prior to each use of the fluorescence application INFRARED 800 with FLOW 800, perform a function test with the INFRARED 800 fluorescence target.
- ▶ Use the fluorescence application INFRARED 800 with FLOW 800 only if this function test has been successfully passed.

Prerequisite

- ☑ KINEVO 900 is switched on and ready for operation
- ☑ The "Fluorescence" function has been assigned [27] to a button on the handgrip or the foot control panel.

Action

1. Place the INFRARED 800 fluorescence target on a flat, level surface.
2. If necessary, remove the attached comparison image with a brief overleaf description from the INFRARED 800 fluorescence target.
3. Focus roughly on the INFRARED 800 fluorescence target.
4. Set the microscope to the desired magnification and correct the image with the fine focusing until it is sharply defined.
5. Check the functioning of the external monitor (option).
6. Actuate the button on the handgrip or foot control panel to which you have assigned the "Fluorescence" function.
 - ⇒ The function test was successful if the image displayed on the touchscreen and on the external monitor corresponds to the attached comparison image of the INFRARED 800 fluorescence target. Slight differences in brightness and inhomogeneities are possible.
 - ⇒ If you activate the "INFRARED 800 Live Overlay" display mode, the fluorescence signal is displayed as a green overlay on the fluorescence target.

5.2 Preparing patients

Prerequisite

- ☑ The tissue to be examined in the surgical field is exposed and can be observed through the KINEVO 900 with good visibility.

Action

- Administer the fluorescent agent approved for the application to the patient only after starting the recording mode.
 - ⇒ The system requires approx. 3 seconds to start the recording mode.

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CAUTION!

Risk of misdiagnosis!

The displayed values are rounded. They serve only to enhance the display and do not represent a measurement function.

- ▶ Do not use the information you obtain from the application INFRARED 800 with FLOW 800 Option for diagnostic purposes.

CAUTION!

Risk of misdiagnosis!

The video monitor is neither calibrated nor designed for diagnostic purposes.

- ▶ Do not use the data displayed on the video monitor for diagnostic purposes.

NOTE

Following failure of the FLOW 800 module, the INFRARED 800 fluorescence application is available.

- ▶ If the FLOW 800 module fails to start following the "Replay" playback phase, you can go on working with the INFRARED 800 fluorescence application as usual.

NOTE

Incorrect presentation of the FLOW 800 display!

Tissue manipulation or setting adjustments on the KINEVO 900 (zoom, focus, position) can give rise to artifacts in the FLOW 800 display. These types of displays are then no longer suitable for use as aids in interpreting the blood flow.

- ▶ Avoid inserting instruments into the operation field during the INFRARED 800 video recording.
- ▶ Do not make any settings adjustments on the KINEVO 900 (zoom, focus, position) during the INFRARED 800 video recording.
- ▶ Do not move the patient while the INFRARED 800 video recording is in progress.
- ▶ If needed, use the "Optimize View" function to remove artifacts.

NOTE**Shaky FLOW 800 display with decrease in brightness at the edges of the image.**

Slight vibrations or uneven illumination of the surgical field during the INFRARED 800 video recording can be corrected by the software.

- ▶ Tap the [Auto-correction] button in the "FLOW 800" menu.

NOTE**Malfunction of the FLOW 800 application due to operating errors!**

The FLOW 800 display will contain artifacts and will not be able to be reasonably interpreted if you change setting parameters and/or move the KINEVO 900 during the INFRARED 800 video recording.

- ▶ Do not change setting parameters such as the focus, zoom, illumination or light field diaphragm during the INFRARED 800 video recording.
- ▶ Do not release the brakes or move the INFRARED 800 during the KINEVO 900 video recording.

You have two options for starting and operating the INFRARED 800 with FLOW 800 Option application. You can combine these two options during the operating sequence:

Action

- ▶ Operating INFRARED 800 with FLOW 800 Option via handgrip or foot control panel [37].
 - ⇒ Actuate the preconfigured buttons on the handgrip or on the foot control panel.
 - ⇒ Briefly press the button on the handgrip or the foot control panel: The application will jump to the next phase.
 - ⇒ Press and hold the button on the handgrip or the foot control panel (longer than 3 seconds): The INFRARED 800 with FLOW 800 Option application will stop and be aborted. In this case, the recorded video data are not saved.
- ▶ Operating INFRARED 800 with FLOW 800 Option via the user interface on the monitor [39].
 - ⇒ Tap the [Fluorescence]  button on the screen and select the fluorescence mode.
 - ⇒ Tap the [Start Record] button: The INFRARED 800 with FLOW 800 Option application will start.
 - ⇒ Switch to the next phase: Tap on the fields at the top of the screen or tap the [Stop Record] button.
 - ⇒ Tap the [Home] button: The INFRARED 800 with FLOW 800 Option application will stop and be aborted. In this case, the FLOW 800 calculations are not generated.

6.1 Operating INFRARED 800 with FLOW 800 Option via handgrip or foot control panel

You can end the fluorescence application at any time by pressing the preconfigured [Fluorescence] button on the hand grip or foot control panel, and holding it down for longer than 3 seconds. INFRARED 800 and FLOW 800 data are lost in the process.

Prerequisite

- ☑ The KINEVO 900 is switched on and ready for operation.
- ☑ The "Fluorescence" function has been assigned [27] to a button on the handgrip or the foot control panel
- ☑ The function test was successfully completed [33]
- ☑ The patient is prepared for the fluorescence application [34]

Action

1. If you have activated the "Skip Setup" function, briefly press the preconfigured [Fluorescence] button on the handgrip or the foot control panel once.
 - ⇒ The "Setup" phase is skipped, the "Record" phase is started and the video recording begins.
2. If you have deactivated the "Skip Setup" function, briefly press the preconfigured [Fluorescence] button on the handgrip or the foot control panel once.
 - ⇒ The screen changes to the setup phase of the INFRARED 800.
 - ⇒ The sequence menu with the different phases of the fluorescence application appears at the top of the INFRARED 800 screen. The blue bar underneath the sequence menu marks the current phase.
 - ⇒ Two broken vertical lines appear which mark the area subsequently displayed in the FLOW 800 module.
 - ⇒ Information indicators for the "Focus" and "Zoom" functions appear on the left edge of the screen.
3. Check the coloring of the two "Focus" and "Zoom" information indicators.
 - ⇒ Red indicator lines indicate that the values set are outside of the recommended range. The values can be adjusted via the corresponding buttons on the handgrip or foot control panel until they are within the recommended range (white indicator lines).
4. If the additional illumination is switched on on the KINEVO 900, a message will appear underneath the two slide controls stating that additional illumination disturbs the fluorescence application.

NOTE! In order for the FLOW 800 display to be illuminated as evenly as possible, the additional illumination should be switched off!

To do this, exit the setup phase by pressing the [Fluorescence] button once on the handgrip or foot control panel. Hold the button longer than 3 seconds to end the INFRARED 800 sequence. Deactivate the additional illumination in the main menu bar [16] of the KINEVO 900 user interface. Restart the INFRARED 800 sequence.

5. Briefly press the [Fluorescence] button on the handgrip or foot control panel again before the fluorescent agent is administered.

⇒ The "Record" phase will start and the video recording will begin.

6. In the recording phase, you can readjust the gain regardless of whether the "AutoGain" function is activated or deactivated by pressing the two yellow buttons on the right handgrip.

7. To end the recording phase, briefly press the [Fluorescence] button on the handgrip or foot control panel again.
The recording duration is limited to a maximum of 5 minutes, i.e. the device automatically ends the recording after 5 minutes.

⇒ The "Record" phase is terminated.

⇒ The "Replay" phase begins: First the short replays are played back. Then the long replays are played back.

8. Navigate through the replays using the joystick on the right handgrip:

Press the joystick up: Playback

Press the joystick down: Pause

Press the joystick to the side: Jump forward/backward one video frame

9. The "Replay" phase is terminated when all short and long replays have been played back. However, if you would like to abort the "Replay" phase, briefly press the [Fluorescence] button on the handgrip or foot control panel again.

⇒ If you have a license for the INFRARED 800 with FLOW 800 application, the FLOW 800 menu will open automatically as soon as the displays have been calculated.

⇒ If you have a license for the INFRARED 800 or INFRARED 800 Compact application, the fluorescence application ends at this point.

Note:

If you have a license for the "PositionMemory" application, the system will automatically save this position. With the following fluorescence applications, you can return to this position and make another INFRARED 800 video recording with the same settings and

the same OPMI position. Detailed information on the "Position-Memory" function is provided in the instructions for use of the KINEVO 900.

This supports subsequent comparisons, especially for FLOW 800.

6.2 Operating INFRARED 800 with FLOW 800 Option via the monitor

You can end the fluorescence application at any time by pressing the preconfigured [Fluorescence] button on the hand grip or foot control panel once and holding it down for longer than 3 seconds or by tapping on the  button on the main screen.


Figure 7: Starting the INFRARED 800 with FLOW 800 Option application

1	[Fluorescence] button	2	[INFRARED 800 with FLOW 800] button
3	[FLOW 800] button	4	[Start Record] button

Prerequisite

- ☒ The KINEVO 900 is switched on and ready for operation.
- ☒ The "Fluorescence" function has been assigned [27] to a button on the handgrip or the foot control panel
- ☒ The function test was successfully completed [33]
- ☒ The patient is prepared for the fluorescence application [34]

Action

1. Tap on the [Fluorescence]  button on the monitor.
2. Select the desired fluorescence mode (INFRARED 800 or INFRARED 800 with FLOW 800).
3. If you have activated the "Skip Setup" function, tap the [Start Record] button.
 - ⇒ The "Setup" phase is skipped, the "Record" phase is started and the video recording begins.

4. If you have deactivated the "Skip Setup" function, tap the [Start Record] button.
 - ⇒ The screen changes to the setup phase of the INFRARED 800.
 - ⇒ The sequence menu with the different phases of the fluorescence application appears at the top of the INFRARED 800 screen. The blue bar underneath the sequence menu marks the current phase.
 - ⇒ Two broken vertical lines appear which mark the area subsequently displayed in the FLOW 800 module.
 - ⇒ The [IR 800 Video] field with the [Start / Stop Record] button appears at the bottom of the screen.
 - ⇒ Information indicators for the recommended "Focus" and "Zoom" settings appear on the left edge of the screen.
5. Check the coloring of the two "Focus" and "Zoom" information indicators.
 - ⇒ Red indicator lines indicate that the values set are outside of the recommended range. The values can be adjusted via the corresponding buttons on the handgrip or foot control panel until they are within the recommended range (white indicator lines).
6. If the additional illumination is switched on on the KINEVO 900, a message will appear underneath the two slide controls stating that additional illumination disturbs the fluorescence application.

NOTE! In order for the FLOW 800 display to be illuminated as evenly as possible, the additional illumination should be switched off!

To do this, exit the setup phase by pressing the [Fluorescence] button once on the handgrip or foot control panel. Hold the button longer than 3 seconds to end the INFRARED 800 sequence. Deactivate the additional illumination in the main menu bar [16] of the KINEVO 900 user interface. Restart the INFRARED 800 sequence.
7. Tap the [Start Record] button in the [IR 800 Video] field before the fluorescent agent is administered.
 - ⇒ The "Record" phase will start and the video recording will begin.
8. To end the recording phase, tap the [Short Replay] button in the sequence menu at the top of the screen or tap the [Stop Record] button in the [IR 800 Video] field.

The recording duration is limited to a maximum of 5 minutes, i.e. the device automatically ends the recording after 5 minutes.

 - ⇒ The "Record" phase is terminated.
 - ⇒ The "Replay" phase begins: First the short replays are played back. Then the long replays are played back.

9. Use the [Start/Stop], [Fast-forward] and [Rewind] control buttons on the monitor to navigate through the replays.
10. The "Replay" phase is terminated when all short and long replays have been played back.
As soon as the FLOW 800 data have been calculated, you can abort the "Replay" phase by tapping the [FLOW 800] button at the top of the screen.
 - ⇒ If you have a license for the INFRARED 800 with FLOW 800 application, the FLOW 800 menu will open automatically as soon as the displays have been calculated.
 - ⇒ If you only have a license for the INFRARED 800 or INFRARED 800 Compact application, the fluorescence application ends at this point.

Note:

If you have a license for the "PositionMemory" application, the system will automatically save this position. With the following fluorescence applications, you can return to this position and make another INFRARED 800 video recording with the same settings and the same OPMI position. Detailed information on the "Position-Memory" function is provided in the instructions for use of the KINEVO 900.

This supports subsequent comparisons, especially for FLOW 800.

6.3 Analyzing and editing FLOW 800 data

6.3.1 Calling FLOW 800 data

There are two ways for you to open the video data of your patient recorded by the INFRARED 800 fluorescence application and processed by the FLOW 800 module:

- Immediately after the INFRARED 800 fluorescence application and all FLOW 800 displays have been calculated [43]
- Via the "Patient" menu [42]

6.3.1.1 Calling up FLOW 800 data via the "Patient" menu

Action

1. Tap on the [Patient] button in the status bar.
⇒ The patient list is displayed.
2. Tap on the patient whose FLOW 800 data you would like to view and/or edit.
⇒ This patient is displayed as the active patient (highlighted in blue) in the status bar.
3. The FLOW 800 data are saved in the patient image data under the "Data" sort key. Tap on the [Data] button in the preview window of the patient image data.

- ⇒ All FLOW 800 data of the selected patient are displayed in the preview window of the patient image data.
4. Tap on the desired preview image.
⇒ The FLOW 800 data are loaded.
⇒ Now you can observe and analyze the video data of the fluorescing tissue precisely.

6.3.1.2 Viewing the FLOW 800 display after the INFRARED 800 fluorescence application

Prerequisite

- ☑ The INFRARED 800 fluorescence application has been completed, the FLOW 800 display has been calculated.

Action

- The FLOW 800 display will automatically appear on the screen after the "Replay" playback phase.
 - ⇒ The FLOW 800 display has been calculated and loaded.
 - ⇒ Now you can observe and analyze the INFRARED 800 video prepared with FLOW 800.

6.3.2 Observing different characteristics of the FLOW 800 display


Figure 8: Characteristics of the FLOW 800 display

1	[Delay] button	2	[Intensity] button
3	[Speed] button	4	[Photo of current display] button
5	Starting point of "color coding" slide control	6	End point of "color coding" slide control

Prerequisite

- ☑ The INFRARED 800 fluorescence application has been completed, the FLOW 800 display has been calculated
- ☑ The "FLOW 800" menu has been opened [42]

Action

1. The "FLOW 800" menu starts with the display of the "Delay" characteristic.

- ⇒ This characteristic displays the time course (in sec.) of the fluorescence influx into the tissue. The red colored areas are perfused first; the blue colored areas are perfused later (delayed).
- 2. Tap the [Intensity] button or press the joystick on the right handgrip once (to the right) to display the "Intensity" characteristic.
 - ⇒ This characteristic displays the maximum fluorescence achieved in the tissue. Values are shown as % of maximum value (100%). The fluorescence intensity is high in the red colored areas and low in the blue colored areas.
- 3. Tap the [Speed] button or press the joystick on the right handgrip to display the "Speed" characteristic.
 - ⇒ This characteristic displays the speed (AI/sec.) of the brightness rise in the tissue. The rise rate is high in red colored areas and low in blue colored areas.
- 4. A slide control for individually adjusting the color coding is assigned to each characteristic on the left edge of the screen. You can randomly and continuously adjust the starting and end points of the slide controls in order to individually adapt the color coding for each characteristic.
 - ⇒ Adjustment of the slide controls changes the presentation of the FLOW 800 display.
- 5. You can generate and save documentation for each characteristic and each setting in the current display. Tap the [Photo]  button at the bottom right of the screen or the correspondingly configured button on the handgrip or foot control panel.
 - ⇒ The documentation of the current display is saved to the patient image data. You can open and view these patient image data via the "Patient" menu at any time.

6.3.3 Setting markers (ROIs)

In the "FLOW 800" menu, you can set up to eight markers (ROIs or regions of interest) at locations in the tissue which are of interest to you.

The diagram displays the time course of the fluorescence intensity for the previously set markers (ROIs). You can observe the change of intensity of the fluorescence brightness (in %) over time (in sec.) in the diagram. Each marker (ROI) is displayed as a separate curve in the diagram [47].


Figure 9: Setting and editing markers (ROIs)

1	Active marker	2	[Show/hide ROI] button
3	[Delete] button	4	"Options" submenu
5	[Auto-correction] button	6	[Photo of current display] button

Prerequisite

- ☑ The INFRARED 800 fluorescence application has been completed, the FLOW 800 display has been calculated.
- ☑ The "FLOW 800" menu has been opened [42]

Action

1. Tap on the location in the FLOW 800 display of interest to you.
 - ⇒ A colored outline (square or adaptive form) then appears at the location you tapped on, the associated value at this marked location is shown underneath the outline and highlighted in blue.
2. You can change the position of each marker at any time by tapping on the respective colored outline and sliding it on the screen.

3. Open the additional "ROI" menu by tapping on the arrow key located in the center on the right edge of the screen.
 - ⇒ The additional "ROI" menu opens on the right-hand side of the screen.
 - ⇒ The currently active marker is highlighted in blue.
4. To hide a marker, tap on the respective marker on the FLOW 800 display or in the additional "ROI" menu.
Then tap the [Show/hide ROI]  button.
 - ⇒ The marker is hidden, but not deleted.
 - ⇒ The marker is displayed "crossed out": 
5. To show a hidden marker again, tap on the "crossed-out" marker in the additional "ROI" menu.
Then tap the [Show/hide ROI]  button.
 - ⇒ The marker is shown again on the FLOW 800 display.
 - ⇒ The marker is displayed as a colored square in the additional "ROI" menu.
6. To delete a marker, tap on the respective marker on the FLOW 800 display or in the additional "ROI" menu.
Then tap on the [Delete]  button.
7. To change the color of a marker, open the "Options" submenu. Tap on the respective marker on the FLOW 800 display or in the additional "ROI" menu. Now select the desired appearance for the marker:
 - Size
 - Shape (square or adaptive)
 - Color
8. Tap the [Auto-correction] button to compensate for and correct small movements of the KINEVO 900 (e.g. vibrations) and, if necessary, poor illumination towards the edges of the image (vignetting).
9. If you have already prepared one or more FLOW 800 displays of your patient and have already set markers there, you can transfer [54] these markers to your current FLOW 800 display.
10. You can generate and save documentation for each diagram and each diagram setting in the current display. Tap the [Photo]  button at the bottom right of the screen or the correspondingly configured button on the handgrip or foot control panel.
 - ⇒ The documentation of the current display is saved to the patient image data. You can open and view these patient image data via the "Patient" menu at any time.

6.3.4 Interpreting and editing the diagram

The diagram displays the time course of the fluorescence intensity for the previously set markers (ROIs). You can observe the change of intensity of the fluorescence brightness (in %) over time (in sec.) in the diagram. Each marker (ROI) is displayed as a separate curve in the diagram.

In addition, the following values are displayed:

- "Delay" characteristic
- "Speed of fluorescence increase" characteristic
- Time to Peak (TTP) (Period of time the contrast medium requires to reach its maximum concentration in the tissue (time until maximum hyperintensity))
- Rise Time (RT) (Period of time the contrast medium requires until an increase from a preset start value to a preset maximum value occurs)


Figure 10: Evaluating a diagram

1	Current playback position	2	Display of parameters
3	Additional "ROI" menu	4	Video recording
5	Video recording slide control	6	[Fast-forward] / [Rewind] button
7	[Start] / [Pause] button	8	[Photo of current display] button
9	"Time" end point slide control	10	"Time" starting point slide control
11	"Intensity" starting point slide control	12	"Intensity" end point slide control

Prerequisite

- ☑ The INFRARED 800 fluorescence application has been completed, the FLOW 800 data have been calculated.
- ☑ The "FLOW 800" menu has been opened [42]
- ☑ At least one marker (ROI) was set [45]

Action

1. Tap the [Diagram] button in the "FLOW 800" menu.
 - ⇒ The "Diagram" menu opens.
 - ⇒ A diagram displaying the markers you set as curves appears on the screen.
 - ⇒ The corresponding INFRARED 800 video is displayed to the right of the diagram.
2. Tap the [Start]  button to play back the video recording. The video recording runs repeatedly and successively until you tap the [Pause]  button.
 - ⇒ The point of time currently displayed in the video is indicated by a white bar in the diagram. This enables you to compare the video display with the diagram display.
3. Tap the [Fast-forward]  and [Rewind]  buttons to move the video recording one sequence forward or back. You can also move the slide control underneath the video recording to jump to a specific location.
 - ⇒ A vertical white bar moves across the diagram, thus marking the current playback position of the video recording.
4. You can adjust the axes of the diagram individually by shifting the starting and end points of the two slide controls (to the left of and underneath the diagram).
 - ⇒ The section from the diagram can be enlarged and displayed in detail by minimizing the areas for "Intensity" (vertical slide control) and "Time" (horizontal slide control).
 - ⇒ The minimized section is marked by two small, white bars on the slide control underneath the video recording. If you now play back the video recording, only this section will be displayed.
5. In order to display the entire diagram, you must reset the starting and end points of both slide controls.
 - ⇒ The entire diagram is then displayed.
 - ⇒ The complete video recording can once again be played back.
6. The "Intensity" characteristic is displayed in the diagram for each marker. If you would like to show all parameters, tap on the desired marker in the additional "ROI" menu.
 - ⇒ All other parameters are shown to the right of the diagram.

- ⇒ The respective curve is highlighted in the diagram and the frame surrounding the parameters has the same color as the selected marker.
- 7. You can generate and save documentation for each diagram and each diagram setting in the current display. Tap the [Photo]  button at the bottom right of the screen or the correspondingly configured button on the handgrip or foot control panel.
 - ⇒ The documentation of the current display is saved to the patient image data. You can open and view these patient image data via the "Patient" menu at any time.
 - ⇒ If you have taken several photos, you can compare these photos (however, only photos of one and the same patient) [52].

6.3.5 Optimize View

Using this function, you can select a partial sequence of the INFRARED 800 video recording that is to be summarized in the FLOW 800 module. In this way you can e.g. exclude artifacts which occurred at the beginning or at the end of the INFRARED 800 video recording from the summary, thus optimizing the display. You then can save the optimized display; the entire, unedited video recording will nevertheless be preserved.


Figure 11: Structure of the "Optimize View" menu

1	Video recording	2	Diagram
3	[Fast-forward] / [Rewind] button	4	[Cut out] button
5	[Cut in] button	6	Timeline selection window end bar
7	Time marker in timeline selection window	8	Video recording slide control
9	Timeline selection window start bar	10	[Start] / [Stop] button

Prerequisite

- ☒ The INFRARED 800 fluorescence application has been completed, the FLOW 800 display has been calculated.
- ☒ The "FLOW 800" menu has been opened [42]
- ☒ At least one marker (ROI) was set [45]

Action

1. Tap the [Optimize View] button.
⇒ A new menu for selecting a partial sequence opens.
2. Tap the [Start]  button to play back the video recording.

3. Tap the [Fast-forward]  and [Rewind]  buttons to move the video recording one sequence forward or back. You can also move the slide control underneath the video recording to jump to a specific location.
4. If you see a disturbance in the beginning of the INFRARED 800 video recording which you would like to exclude from the summary, tap the [Pause] button  as soon as this disturbance has been played back and the normal fluorescence recording is visible.
 - ⇒ The current playback position in the video recording is displayed in the timeline selection window as a time marker (small white bar).
5. Drag the blue start bar of the timeline selection window onto the time marker.
6. Tap the [Cut in] button.
 - ⇒ The FLOW 800 display starts at this point.
 - ⇒ In the diagram view to the right of the video recording, the cut-out area is colored black. The dark blue colored area marks the remaining playback time.
7. Tap the [Start]  button to continue the playback of the video recording.
8. Tap the [Pause] button  if you do not want to play back the entire video recording, but would like to remove the later portion of the video from the FLOW 800 display.
 - ⇒ The current playback position in the video recording is displayed as a time marker (small white bar) in the timeline selection window.
9. Drag the blue end bar of the timeline selection window onto the time marker.
10. Tap the [Cut out] button.
 - ⇒ The FLOW 800 display ends at this point of time.
 - ⇒ In the diagram view to the right of the video recording, the cut-out area is colored black. The dark blue colored area marks the remaining playback time.
11. You can move the start and end points of the INFRARED 800 video used for the FLOW 800 display at any time as long as you have not tapped the [Apply] button.
12. To apply the optimized display, tap the [Apply] button.
 - ⇒ The selected partial sequence will be used to generate an optimized display of the FLOW 800 data.

- ⇒ The optimized display is used only temporarily. You can take a photo of the optimized display. Tap the [Photo]  button at the bottom right of the screen or the correspondingly configured button on the handgrip or foot control panel.
- ⇒ The original FLOW 800 displays remain unaffected and can be restored at any time.
- ⇒ The "Optimize View" menu is closed.
- ⇒ The "Diagram" menu opens and the optimized, edited display is shown.

13. You can play back the abridged video recording.

14. Tap the [Optimize View] button to edit the display again. Adapt your display gain if necessary.

6.3.6 Comparing FLOW 800 displays

You can compare the diagram view of the current recording with the diagram views of previously created recordings of the same patient, provided that you previously saved the diagram views of the other recordings.

For example, you can record INFRARED 800 videos of your patient at the beginning and at the end of the surgical procedure in order to directly compare the fluorescence intensity in the tissue based on the Delay, Intensity and Speed characteristics at the beginning and at the end of surgery.

You can compare different diagrams only if you have generated and saved [47] a documentation of the current display for each previously created diagram.

Prerequisite

- ☒ The INFRARED 800 fluorescence application has been completed, the FLOW 800 display has been calculated.
- ☒ At least two FLOW 800 displays of the same patient are available.
- ☒ The "FLOW 800" menu is open

Action

1. Start the "FLOW 800" [42] menu to compare your currently acquired FLOW 800 display with a display acquired at another time or open a FLOW 800 data set in the "Patient" menu.
2. Tap the [Compare View] slide switch.
 - ⇒ The slide switch lights up in blue .
 - ⇒ The "Compare View" menu opens.
 - ⇒ The current FLOW 800 display is displayed on the right half of the screen.

- ⇒ All available FLOW 800 displays for a given patient are displayed as preview images at the bottom of the left half of the screen.
- 3. Select a FLOW 800 display from the preview images and tap on it.
 - ⇒ The selected FLOW 800 diagram is displayed in the left half of the screen.
 - ⇒ The preview image of the selected FLOW 800 display is highlighted in blue.
- 4. You can view the different characteristics by tapping the corresponding buttons, i.e. [Delay], [Intensity] and [Speed].
 - ⇒ The images of previously and currently recorded FLOW 800 displays are displayed side by side.
- 5. If you have already documented diagrams for the current patient, these can be compared with the current diagram. To do this, tap the [Diagram] button.
 - ⇒ The current diagram is displayed on the right half of the screen.
 - ⇒ The saved diagrams are displayed as preview images in the bottom left half of the screen.
- 6. Select a diagram from the preview images and tap on it.
 - ⇒ The selected diagram is displayed in the left half of the screen.
 - ⇒ The preview image of the selected diagram is highlighted in blue.
- 7. Deactivate compare mode by tapping on the [Compare View] slide switch again.

6.3.7 Transferring markers (ROIs)

This function allows you to apply the markers from a INFRARED 800 video that has already been analyzed with FLOW 800.

Prerequisite

- ☑ A previously recorded INFRARED 800 video was analyzed with FLOW 800, and markers (ROIs) were set

Action

1. Start a new INFRARED 800 video recording.
 - ⇒ Once the INFRARED 800 video recording has ended, the "FLOW 800" menu will open.
 - ⇒ The "FLOW 800" menu on the right side of the screen contains the [Apply ROIs from earlier analysis] button.
2. Tap the [Apply ROIs from earlier analysis] button.
 - ⇒ The markers from the previous FLOW 800 data set are applied to the current FLOW 800 display.
3. Reposition the ROIs in your current FLOW 800 display, if needed.
 - ⇒ The size, color and positions of the markers are applied and are readily available for analysis purposes.
4. You can view the different characteristics by tapping the corresponding buttons, i.e. [Delay], [Intensity] and [Speed].
5. Tap the [Compare View] slide switch if you would like to view or compare the markers of the two FLOW 800 displays side by side.
 - ⇒ The two FLOW 800 displays are shown side by side on the screen. This enables you to compare the markers directly.
6. To quit the "FLOW 800" menu, tap the  button.

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7.5 Ambient requirements for operation

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7.6 Ambient requirements for transport and storage

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7.7 Operation with intraoperative MRT systems

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Glossary

AutoDetection

The system determines the point of time in the video data when the fluorescence agent flows into the tissue. If this function is switched on, the playback phase automatically starts at this point of time. Video data which precede this point of time are not displayed.

Automatic Gain Control

Automatic camera gain increase during the recording phase.

AVM

Arteriovenous malformation (vascular malformations of the brain)

Rise Time (RT)

Period of time in which the fluorescence intensity increases from 10% to 90% of the maximum value achieved.

ROI (region of interest)

During the playback of the video data, you can set markers (ROIs: regions of interest). The intensity of the fluorescence in the tissue is processed and displayed as a curve in a diagram for each marker set. With the help of this diagram, you can conveniently and quickly evaluate and compare the time course of the fluorescence intensity in the tissue.

Time to Peak (TTP)

Period of time the contrast medium requires to reach its maximum concentration in the tissue (time until maximum hyperintensity)

Empty page, for your notes

Keyword index

A

Access the help menu	16
Adjusting the screen display	28
Auto-correction	
Description	14
Starting	46
AutoDetection	
Description	13
Powering the system up/down	29
AutoGain	
Description	13
Powering the system up/down	29

C

Characteristics	
Description	13
Observing and interpreting	43
Configuring the [Fluorescence] button	27

D

Date, time (display)	16
Diagram	
Interpreting and editing	47
Menu description	22
Documenting current display	43

F

FLOW 800	
Applying markers	54
Characteristics (description)	13
Comparing displays	52
Documenting current display	43
Functional description	12
Hand grip/foot control panel operation ..	37
Menu description	20

Observing characteristics	43
Opening data	42
Operation via monitor	39
Setting markers	45
Showing video data in eyepieces	31

Fluorescence

Button assignment	27
Setting the parameters	28

I

INFRARED 800

Functional description	12
Menu description	18
Record (recording phase)	13
Replay (playback phase)	13
Setup (preparation phase)	13
Showing video data in eyepieces	31

INFRARED 800

Hand grip/foot control panel operation ..	37
Operation via monitor	39

M

Main menu bar	17
Markers	
Applying	54
Setting	45

O

Optimize View	
Menu description	24
Optimizing display	
Editing	50

R

Replay (playback phase)	
-------------------------	--

Long replay (description)	13
Setting a long replay	30
Setting a short replay	30
Short replay (description)	13

S

Status bar	16
------------------	----

T

Target group	9
--------------------	---

U**User Interface**

"Diagram" menu	22
"Optimize View" menu	24
"FLOW 800" menu	20
"INFRARED 800" menu	18
Access the "Extended operation" menu	16
Access the configuration menu	16
Main menu description	16
Open the "Help" menu	16

V**Video recording**

Editing	50
Start/stop via hand grip or foot control panel	37
Start/stop via monitor	39

Empty page, for your notes

Carl Zeiss Meditec AG

Goeschwitzer Strasse 51-52

07745 Jena

Germany

www.meditec.zeiss.com

CE 0297

Subject to change.

G-30-1956-en - 4.0 - 2020-03-19

FOR THE RUSSIAN FEDERATION ONLY

The instruction for use «Fluorescent module INFRARED 800 with FLOW 800 option. For KINEVO 900» (Carl Zeiss Meditec AG) in accordance with national and international requirements, is standard for all countries. In accordance with national requirements and in order to comply with Russian legislation the Instruction for use is supplemented by these Addendums, which are the parts of the Instruction for use for the products supplied to the Russian Federation.

Addendum I to the document Instruction for use «Fluorescent module INFRARED 800 with FLOW 800 option. For KINEVO 900»

Paragraph 2.2.2 «Normal use» (page 19) should be supplemented with the following information:

The fluorescent substances required for operation are not included in the scope of product delivery, but are supplied separately by other companies. Responsibility for quality, efficiency and safety of the fluorescent substances, state registration and compatibility with fluorescence module INFRARED with FLOW 800 option lies on the manufacturer of the corresponding substances.

Only use the fluorescent substances that are registered in RF.

Paragraph 2 «Safety notes» should be supplemented with sub-paragraph 2.3 «Contraindications»:

INFRARED 800 with FLOW 800 fluorescent module
provided that ZEISS module is used within its indication of use, there are no known contraindications.

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