IFU Vacuum Tubes
Sodium Citrate

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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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
Signature:

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(Date: Jul. 21, 2020)

证明书验证网站: Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

INSTRUCTION FOR USE

Signature, Stamp

Approved
by: BAI QIANG
Position: SALES MANAGER

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Vacuum tube VacPlus with Sodium Citrate 3,2% (blue)
Vacuum tube VacPlus with Sodium Citrate 3,2% (white)
Vacuum tube VacPlus with Sodium Citrate 3,8% (blue)
Vacuum tube VacPlus with Sodium Citrate 3,8% ESR (black)
RC № XXX от XX.XX.XXXX

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This product consists of a plastic or glass¹ tube, a rubber stopper, a protective cap, a label and additive: Sodium Citrate 3,2%, 3,8% and Sodium Citrate 3,8% ESR.

¹Vacuum tube VacPlus with Sodium Citrate 3,8% ESR (black) sizes 8x120 9x120 is made of glass. Other sizes of this tube are made of PET.

is made of PET, which practically does not break and is a clinically tested, safer alternative to the stopper tightly closes the tube, regardless of the presence of a vacuum in it. The protective cap is designed as a plug cap for easy connection and detachment, which eliminates the possibility of contact with the health worker's hands and blood samples that may cause sample contamination.

Test tubes are sterile, equipped with a convenient label with a field for recording with double tear-off digital code. Each tube has a mark to accurately determine the level of filling (filling line) located on the label or on the tube itself, depending on the intended volume. Compliance with the required level of filling allows to achieve the correct ratio of blood and filler (additive). Test tubes are available in a wide range of volumes. Test tube sizes are suitable for most centrifuge models. The inner walls of the tubes are covered with fillers (additives) applied by micro dispersed spraying. The caps of the test tubes are color-coded in accordance with the international standard ISO 6710.

Table 1

Name of the tube	Cap (plug) color
Vacuum tube VacPlus with Sodium Citrate 3,2%	Blue
Vacuum tube VacPlus with Sodium Citrate 3,2%	White
Vacuum tube VacPlus with Sodium Citrate 3,8%	Blue
Vacuum tube VacPlus with Sodium Citrate 3,8% ESR	Black

Blood coagulation tubes contain a buffered solution of trisodium citrate. Citrate concentrations of 0.109 mol / L (3.2%) or 0.129 mol / L (3.8%) are available. The choice of concentration depends on the laboratory policy. The blood reagent ratio is 9: 1.

ESR tubes contain a 3.8% buffered solution of trisodium citrate (0.129 mol / L). The ratio of blood reagent is 4: 1. Tubes for determining ESR are used to collect and transport venous blood and then determine the erythrocyte sedimentation rate.

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Purpose: Test tubes are used for coagulological studies, as well as for the study of platelet function (research platelet-rich plasma) and coagulation factors for deep plasma freezing

Test tubes for ESR are used to determine the erythrocyte sedimentation rate (ESR) well as for the study of platelet function (research platelet-rich plasma)

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PRECAUTIONS. WARNINGS. CAUTIONS

Precautions

1. Handle all biological samples and sharp objects in accordance with the rules adopted at your medical facility
2. Discard used sharps into a special container for disposal.
3. All liquid reagents are clear. Do not use tubes if the additives have changed color or there is a foreign substance in them.

not use tubes after the expiration date.

qs:

- Please check the directions for use and designation.
- Check the tube for damage, contamination, leakage.
- Verify that the tube specification and test parameter match.

Precautions and storage requirements (indicated on the label):

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INSTRUCTION FOR USE

- 1) Identify the patient, introduce himself, and explain the course and purpose of the procedure. Ensure that the patient has informed consent for the upcoming blood collection procedure. In the absence of such, clarify further actions with a doctor.
- 2) Offer the patient or help him to take a comfortable position: sitting or lying down.
- 3) Mark the tubes, indicating the patient's full name, department (in order to eliminate errors in the identification of the biomaterial sample). Using the digital code make appropriate notes in the logbook. (Be sure to put on the label on the tube, the name of the medical organization, the name and initials of the patient, his gender, age, date and time of blood collection and send the tube to the laboratory. Labels should always be filled in the presence of the patient.)
- 4) Treat your hands hygienically, dry them.
- 5) Treat your hands with an antiseptic. Do not dry, wait until the antiseptic is completely dry.
- 6) Wear medical gloves.
- 7) Select tubes corresponding to the declared tests or laboratory tests required by the patient, prepare a needle, holder, alcohol wipes, and a plaster.
- 8) Place a tourniquet on a shirt or diaper 7-10 cm above the venipuncture site. Harness need to impose no more than one minute. A longer vein squeeze time may affect the test results due to changes in the concentration of certain parameters in the blood.
- 9) Ask the patient to make a fist. It is impossible to set for a hand an exercise stress (vigorous "squeezing and unclenching a fist"), since this may lead to changes in blood concentrations of some indicators. Choose venipuncture site. The most common are the ulnar and saphenous veins, but smaller and fuller veins of the dorsal surface of the wrist.

Performing the procedure:

- 1) Take the needle and remove the protective cap from it. If a double-sided needle is used, remove the protective cap in gray or white.
- 2) Insert the needle into the needle holder and screw in until it stops.
- 3) Disinfect the venipuncture site with a gauze napkin or swab moistened with an antiseptic solution, in a circular motion from the center to the periphery.
- 4) Wait until the antiseptic solution is completely dry (30-60 s). You cannot wipe and blow the puncture site, so as not to bring microorganisms. It is also impossible to palpate the vein after disinfection. If difficulties arise during venipuncture and the vein is palpated again, this area must be disinfected again.
- 5) Remove the cap on the other side of the needle.
- 6) To clasp the patient's forearm with the left hand so that the thumb is 3-5 cm below the venipuncture site.

in the needle along the line with the vein bevel up and puncture the vein at an angle 15-30 °.

Insert the previously prepared tube into the needle holder until it stops and hold it until the blood flowing into the tube. The harness must be removed immediately after the start of blood flow into the tube. Make sure the patient opens his fist. The blood passes into the tube until it fully compensates for the vacuum created in it. If the blood does not flow, it means that the needle went through the vein - in this case, you need to slightly pull out the needle, but do not remove it until the blood goes to the tube. The accuracy of filling the tube is $\pm 10\%$ of the nominal volume.

9) Remove the tube from the holder.

10) Immediately after filling, the tube must be gently inverted to mix the sample with the filler: a tube without anticoagulants - 5-6 times; test tube with citrate - 3-4 times; a tube with heparin, EDTA and other additives - 8-10 times. The tubes cannot be shaken - it can cause foaming and hemolysis, as well as lead to mechanical lysis of red blood cells.

If necessary, a number of other tubes are inserted into the needle holder to obtain the desired volume of blood for various studies. Re-insert the needle into the vein for this is not necessary.

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on of samples for analysis

Additive	Maximum time after blood sampling	How many times to flip	Recommended relative centrifugal acceleration	Recommended centrifugation time, min
Sodium Citrate 9:1	3 hours	8-10 times	2000-3000 g	10-15
Sodium citrate 4:1 (ESR)	3 hours	6-8 times	-	-

SYMBOLS WITH DESCRIPTION

	Batch code
	Used by date
	Manufacturer
	Date of manufacturer
	Do not reuse
	Do not use if pack is damaged
	Sterilized by radiation
	Keep away from sunlight
	Keep dry
	Temperature limit 4-25°C

	Fragile
	Up
	Layer constraint
	CE marking
	In vitro diagnostic medical device

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Values of the parameters, characteristics (features) of the medical device are presented in the table below. Where applicable, the permissible deviations from the nominal characteristics are $\pm 10\%$, unless otherwise agreed.

Name	Volume	Dimensions	Additive	Cap (plug) color	Quantity
Vacuum tube VacPlus with Sodium Citrate 3,2% (blue). Blood volume: 1.8 ml, 2.7 ml, 3.6 ml, 4.5 ml, 5.4 ml, 6.3 ml, 7.2 ml, 8.1 ml, 9 ml. Dimensions: 13 mm x 75 mm, 13 mm x 100 mm, 16 mm x 100 mm.	1.8 ml, 2.7 ml, 3.6 ml, 4.5 ml, 5.4 ml, 6.3 ml, 7.2 ml, 8.1 ml, 9 ml	13x75 mm, 13x100 mm, 16x100 mm	Sodium citrate 3,2%	Blue	1, 2, 50, 100 pcs/rack
Vacuum tube VacPlus with Sodium Citrate 3,2% (white). Blood volume: 1.8 ml, 2.7 ml, 3.6 ml, 4.5 ml, 5.4 ml, 6.3 ml, 7.2 ml, 8.1 ml, 9 ml. Dimensions: 13 mm x 75 mm, 13 mm x 100 mm, 16 mm x 100 mm.	1.8 ml, 2.7 ml, 3.6 ml, 4.5 ml, 5.4 ml, 6.3 ml, 7.2 ml, 8.1 ml, 9 ml	13x75 mm, 13x100 mm, 16x100 mm	Sodium citrate 3,2%	White	1, 2, 50, 100 pcs/rack
Vacuum tube VacPlus with Sodium Citrate 3,8% (blue). Blood volume: 1.8 ml, 2.7 ml, 3.6 ml, 4.5 ml, 5.4 ml, 6.3 ml, 7.2 ml, 8.1 ml, 9 ml. Dimensions: 13 mm x 75 mm, 13 mm x 100 mm, 16 mm x 100 mm.	1.8 ml, 2.7 ml, 3.6 ml, 4.5 ml, 5.4 ml, 6.3 ml, 7.2 ml, 8.1 ml, 9 ml	13x75 mm, 13x100 mm, 16x100 mm	Sodium citrate 3,8%	Blue	1, 50, 100 pcs/rack

tube VacPlus Sodium Citrate 3,8% (black). Blood Volume: 1.28 ml, 1.6 ml, 1.8 ml, 2 ml, 2.4 ml, 3 ml, 4 ml, 5 ml. Dimensions: 8 mm x 120 mm, 9 mm x 120 mm, 13 mm x 75 mm, 13 mm x 100 mm.	1.28 ml, 1.6ml, 1.8 ml, 2 ml, 2.4 ml, 3 ml, 3.2 ml, 4 ml, 5 ml	8x120 mm, 9x120 mm, 13x75 mm, 13x100 mm, 16x100 mm	Sodium citrate 3,8%	Black	50, 100 pcs/rack
--	--	--	------------------------	-------	---------------------

MAIN RAW MATERIALS

Tube - PET or glass (for some sizes of tubes for ESR)

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Group packaging (racks of 100 or 50 pieces of tubes with sodium citrate) can be additionally packaged in an aluminum thermal bag (in foil) to guarantee the specified shelf life of such tubes, which prevents the accelerated evaporation of citrate.

In addition, tubes with sodium citrate can be packed in individual aluminum packaging of 1 or 2 pcs. (see tab. 2).

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Operation

Temperature	15-30° C
Relative pressure	700 - 1060 hPa
Humidity	35% - 75%

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ubes, including expired ones, should be disposed of as Class A waste and can be disposed of as household waste.

to national regulations when disposing or recycling the product, its components and its packaging.

Food sample is a hidden source of contamination. Therefore, all components of the equipment, consumables and waste substances that come into contact with the sample should be used in accordance with biological safety rules.

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Vacuum tube VacPlus with NaHe
Vacuum tube VacPlus with LiHe
Vacuum tube VacPlus with NaHe and separation gel

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This product consists of a plastic tube, a rubber stopper, a protective cap, a label and additive: NaHe, LiHe, NaHe and gel.

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of the tube	Cap (plug) color
Vacuum tube VacPlus with NaHe	green
Vacuum tube VacPlus with LiHe	
Vacuum tube VacPlus with NaHe and separation gel	

Heparin is an anticoagulant that activates antithrombins and prevents blood coagulation, while not distorting the parameters studied. The blood components in the sample with lithium salt of heparin are stable up to 6 hours.

Caution: The use of these tubes for monitoring drugs and determining the lithium content is prohibited.

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Purpose: Test tubes are used for clinical studies of plasma as well for platelet function test (platelet rich plasma (PRP)).

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INSTRUCTION FOR USE

- 1) Identify the patient, introduce himself, and explain the course and purpose of the procedure. Ensure that the patient has informed consent for the upcoming blood collection procedure. In the absence of such, clarify further actions with a doctor.
- 2) Offer the patient or help him to take a comfortable position: sitting or lying down.
- 3) Mark the tubes, indicating the patient's full name, department (in order to eliminate errors in the identification of the biomaterial sample). Using the digital code make appropriate notes in the logbook.

to put on the label on the tube, the name of the medical organization, the name and initials of the patient, his gender, age, date and time of blood collection and send the tube to the laboratory. Labels should always be filled in the presence of the patient.)

Treat your hands hygienically, dry them.

Treat your hands with an antiseptic. Do not dry, wait until the antiseptic is completely dry.

Wear medical gloves.

7) Select tubes corresponding to the declared tests or laboratory tests required by the patient, prepare a needle, holder, alcohol wipes, and a plaster.

8) Place a tourniquet on a shirt or diaper 7-10 cm above the venipuncture site. Harness need to impose no more than one minute. A longer vein squeeze time may affect the test results due to changes in the concentration of certain parameters in the blood.

9) Ask the patient to make a fist. It is impossible to set for a hand an exercise stress (vigorous "squeezing and unclenching a fist"), since this may lead to changes in blood concentrations of some indicators. Choose venipuncture site. The most common are the ulnar and saphenous veins, but smaller and fuller veins of the dorsal surface of the wrist.

Performing the procedure:

1) Take the needle and remove the protective cap from it. If a double-sided needle is used, remove the protective cap in gray or white.

2) Insert the needle into the needle holder and screw in until it stops.

3) Disinfect the venipuncture site with a gauze napkin or swab moistened with an antiseptic solution, in a circular motion from the center to the periphery.

4) Wait until the antiseptic solution is completely dry (30-60 s). You cannot wipe and blow the puncture site, so as not to bring microorganisms. It is also impossible to palpate the vein after disinfection. If difficulties arise during venipuncture and the vein is palpated again, this area must be disinfected again.

5) Remove the cap on the other side of the needle.

6) To clasp the patient's forearm with the left hand so that the thumb is 3-5 cm below the venipuncture site.

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Make sure that the patient does not have external bleeding in the venipuncture area.

Apply a pressure bandage on the arm or a bactericidal patch.

Disinfect gloved hands with a disinfectant.

Disinfect all used material. Remove gloves, place in a disinfection tank or a waterproof bag /

Container for disposal of class B waste.

Treat your hands hygienically, dry them.

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Additive	Maximum time after blood sampling	How many times to flip	Recommended relative centrifugal acceleration	Recommended centrifugation time, min
LiHe/ NaHe	3 hours	8-10 times	2000 - 3000 g	15
NaHe+gel	3 hours before centrifugation and 48 hours after centrifugation	8-10 times	2200 g	15

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Tubes with sodium heparin and separation gel can be packed in individual aluminum packaging of 1 or 2 pcs. (see tab. 2).

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Transportation of products without blood should be carried out by all types of covered transport in compliance with the conditions and requirements established for this type of transport, at the

temperature from 4 to +25 ° C. During transportation to provide protection from the sun, rain and do not mix with harmful or toxic substances. No special requirement stances.

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THIS IS TO CERTIFY THAT: the seal of HEBEI XINLE SCI & TECH CO., LTD. on the annexed DOCUMENT is genuine.

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Signature:

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
(Date: Jul. 21, 2020)

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Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

INSTRUCTION FOR USE

Approved
Name: BAI QIANG
Position: SALES MANAGER

Signature, Stamp
河北鑫乐医疗器械科技股份有限公司
HEBEI XINLE SCI&TECH CO., LTD

Vacuum tube VacPlus with sodium fluoride and potassium oxalate
Vacuum tube VacPlus with sodium fluoride and EDTA K2
RC № XXX or XX.XX.XXXX

TITLE. DESCRIPTION

1. Vacuum tube VacPlus with sodium fluoride and potassium oxalate (gray). Blood volume: 1 ml, 2 ml, 3 ml, 4 ml, 5 ml, 6 ml, 7 ml, 8 ml, 9 ml. Dimensions: 13 mm x 75 mm, 13 mm x 100 mm, 16 mm x 100 mm.
2. Vacuum tube VacPlus with sodium fluoride and EDTA K2 (gray). Blood volume: 1 ml, 2 ml, 3 ml, 4 ml, 5 ml, 6 ml, 7 ml, 8 ml, 9 ml. Dimensions: 13 mm x 75 mm, 13 mm x 100 mm, 16 mm x 100 mm.

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This product consists of a plastic tube, a rubber stopper, a protective cap, a label and additive: sodium fluoride and potassium oxalate, sodium fluoride and EDTA K2.

The tube is made of PET, which practically does not beat and is a clinically tested, safer alternative to glass. The stopper tightly closes the tube, regardless of the presence of a vacuum in it. The protective cap is designed as a plug cap for easy connection and detachment, which eliminates the possibility of direct contact with the health worker's hands and blood samples that may cause sample contamination.

All tubes are sterile, equipped with a convenient label with a field for recording with double tear-off digital code. Each tube has a mark to accurately determine the level of filling (filling line) located on the label or on the tube itself, depending on the intended volume. Compliance with the required level of filling allows to achieve the correct ratio of blood and filler (additive). Test tubes are available in a wide range of volumes. Test tube sizes are suitable for most centrifuge models. The inner walls of the tubes are covered with fillers (additives) applied by micro dispersed spraying. The caps of the test tubes are color-coded in accordance with the international standard ISO 6710.

Table 1

Name of the tube	Cap (plug) color
------------------	------------------

10 mm tube VacPlus with sodium fluoride and potassium oxalate	gray
10 mm tube VacPlus with sodium fluoride and EDTA K2	

contain anticoagulant potassium oxalate and glucose stabilizer sodium fluoride. The anticoagulant prevents blood coagulation by binding calcium ions, and sodium fluoride inhibits glycolysis by blocking the enzyme enolase. The glucose concentration in the sample is stable up to 6 hours at room temperature. A study of glucose concentration should be carried out as soon as possible, since sodium fluoride, inhibiting one of the last stages of glucose breakdown, is less effective against the initial stages.

Attention! Glucose test tubes should be filled to the volume indicated on them, since excess anticoagulant can cause hemolysis. The contents of the tube after blood collection should be mixed with extreme caution.

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Purpose: Test tubes are used to determine the concentration of glucose, lactate in the blood.

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PRECAUTIONS. WARNINGS. CAUTIONS

Precautions

1. Handle all biological samples and sharp objects in accordance with the rules adopted at your medical facility
2. Discard used sharps into a special container for disposal.
3. All liquid reagents are clear. Do not use tubes if the additives have changed color or there is a foreign substance in them.
4. Do not use tubes after the expiration date.

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- Store at temperature from +4 ° C to 25 ° C
- Do not use if packaging is opened or damaged.
- Do not use after expiration date

INSTRUCTION FOR USE

- Identify the patient, introduce himself, and explain the course and purpose of the procedure. Make sure that the patient has informed consent for the upcoming blood collection procedure. In the absence of such, clarify further actions with a doctor.
- Offer the patient or help him to take a comfortable position: sitting or lying down.
- Mark the tubes, indicating the patient's full name, department (in order to eliminate errors in the identification of the biomaterial sample). Using the digital code make appropriate notes in the logbook. (Be sure to put on the label on the tube, the name of the medical organization, the name and initials of the patient, his gender, age, date and time of blood collection and send the tube to the laboratory. Labels should always be filled in the presence of the patient.)
- 4) Treat your hands hygienically, dry them.
 - 5) Treat your hands with an antiseptic. Do not dry, wait until the antiseptic is completely dry.
 - 6) Wear medical gloves.
 - 7) Select tubes corresponding to the declared tests or laboratory tests required by the patient, prepare a needle, holder, alcohol wipes, and a plaster.
 - 8) Place a tourniquet on a shirt or diaper 7-10 cm above the venipuncture site. Harness need to impose no more than one minute. A longer vein squeeze time may affect the test results due to changes in the concentration of certain parameters in the blood.
 - 9) Ask the patient to make a fist. It is impossible to set for a hand an exercise stress (vigorous "squeezing and unclenching a fist"), since this may lead to changes in blood concentrations of some indicators. Choose venipuncture site. The most common are the ulnar and saphenous veins, but smaller and fuller veins of the dorsal surface of the wrist.

Performing the procedure:

- 1) Take the needle and remove the protective cap from it. If a double-sided needle is used, remove the protective cap in gray or white.
- 2) Insert the needle into the needle holder and screw in until it stops.
- 3) Disinfect the venipuncture site with a gauze napkin or swab moistened with an antiseptic solution, in a circular motion from the center to the periphery.
- 4) Wait until the antiseptic solution is completely dry (30-60 s). You cannot wipe and blow the puncture site, so as not to bring microorganisms. It is also impossible to palpate the vein after disinfection. If difficulties arise during venipuncture and the vein is palpated again, this area must be disinfected again.
- 5) Remove the cap on the other side of the needle.
- 6) To clasp the patient's forearm with the left hand so that the thumb is 3-5 cm below the venipuncture site.
- 7) Position the needle along the line with the vein bevel up and puncture the vein at an angle 15-30 ° to the skin.
- 8) Insert the previously prepared tube into the needle holder until it stops and hold it until the blood stops flowing into the tube. The harness must be removed immediately after the start of blood flow into the tube. Make sure the patient opens his fist. The blood passes into the tube until it fully compensates for the vacuum created in it. If the blood does not flow, it means that the needle went through the vein - in this case, you need to slightly pull out the needle, but do not remove it until the blood goes to the tube. The accuracy of filling the tube is $\pm 10\%$ of the nominal volume.
- 9) Remove the tube from the holder.
- 10) Immediately after filling, the tube must be gently inverted to mix the sample with the filler: a tube without anticoagulants - 5-6 times; test tube with citrate - 3-4 times; a tube with heparin, EDTA and other additives - 8-10 times. The tubes cannot be shaken - it can cause foaming and hemolysis, as well as lead to mechanical lysis of red blood cells.

If necessary, a number of other tubes are inserted into the needle holder to obtain the desired volume for various studies. Re-insert the needle into the vein for this is not necessary.

After all the required tubes are filled, attach a dry sterile cloth to the venipuncture site and remove the needle.

Termination of the procedure:

Place the used needle together with the disposable holder in a sharps container. In order to prevent contact with blood, it is forbidden to disassemble the needle and the holder in your hands!

- 2) Ensure that the patient does not have external bleeding in the venipuncture area.
 - 3) Apply a pressure bandage on the arm or a bactericidal patch.
 - 4) Handle gloved hands with a disinfectant.
 - 5) Disinfect all used material. Remove gloves, place in a disinfection tank or a waterproof bag / container for disposal of class B waste.
 - 6) Treat your hands hygienically, dry them.
 - 7) Refine patient well-being.
 - 8) Make an appropriate record of the results of the service in the medical records or arrange a referral.
 - 9) Arrange the delivery of test tubes with the obtained laboratory material to the laboratory
- Be sure to put on the label on the tube, the name of the medical organization, the patient's name and initials, his gender, age, date and time of blood collection and send the tube to the laboratory. Labels must always be filled in the presence of the patient.

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Centrifuging



If necessary, the centrifugation procedure, make sure that the tubes are inserted into the rotor in such a way that the cap does not rest on the walls of the centrifuge cup, otherwise the cap may come off the tubes - a more stable gel barrier is formed in horizontal centrifuges than in centrifuges with a fixed angle of inclination - high temperatures can have a negative effect on the physical properties of the gel. The most complete separation of serum and plasma occurs at a temperature of 20°C - 22°C.

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Additive	Maximum time after blood sampling	How many times to flip	Recommended relative centrifugal acceleration	Recommended centrifugation time, min
NaF and KOx or EDTA K2	24 hours	8-10 times	1500 g	10

SYMBOLS WITH DESCRIPTION

LOT	Batch code
------------	------------

	Used by date
	Manufacturer
	Date of manufacturer
	Do not reuse
	Do not use if pack is damaged
	Sterilized by radiation
	Keep away from sunlight
	Keep dry
	Temperature limit 4-25°C
	Fragile
	Up
	Layer constraint
	CE marking
	In vitro diagnostic medical device

GUARANTEED BY MANUFACTURER VALUES OF MAIN PARAMETERS, CHARACTERISTICS (FEATURES) OF MEDICAL DEVICE.

of the parameters, characteristics (features) of the medical device are presented in the table. Where applicable, the permissible deviations from the nominal characteristics are $\pm 10\%$, as otherwise agreed.

Name	Volume	Dimensions	Additive	Cap (plug) color	Quantity
Vacuum tube VacPlus with sodium fluoride and potassium oxalate (gray). Blood volume: 1 ml, 2 ml, 3 ml, 4 ml, 5 ml, 6 ml, 7 ml, 8 ml, 9 ml. Dimensions: 13 mm x 75 mm, 13 mm x 100 mm, 16 mm x 100 mm.	1 ml, 2 ml, 3 ml, 4 ml, 5 ml, 6 ml, 7 ml, 8 ml, 9 ml	13x75 mm, 13x100 mm	sodium fluoride and potassium oxalate	Gray	50, 100 pcs/rack
Vacuum tube VacPlus with sodium fluoride and K2EDTA (gray). Blood volume: 1 ml, 2 ml, 3 ml, 4 ml, 5 ml, 6 ml, 7 ml, 8 ml, 9 ml. Dimensions: 13 mm x 75 mm, 13 mm x 100 mm, 16 mm x 100 mm.	1 ml, 2 ml, 3 ml, 4 ml, 5 ml, 6 ml, 7 ml, 8 ml, 9 ml	13x75 mm, 13x100 mm	sodium fluoride and K2EDTA	Gray	50, 100 pcs/rack

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Keep vacuum tubes for taking blood without blood at a temperature of 4-25 ° C in a dry room. Keep away from harmful or poisonous substances.
 Avoid direct sunlight. Exceeding the recommended storage temperature may lead to a loss of quality of the mini tube (for example, evaporation or staining of fillers (additives)).

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Shelf life is 24 months

WARRANTY OBLIGATION

All questions related to the operation and quality of the product can be addressed to the Authorized representative of the manufacturer in the territory of the Russian Federation.

The manufacturer guarantees that the product meets the requirements of the standards and documentation of the manufacturer under the conditions of transportation, storage, operation.

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REAL CARE LLC

Perovskayast. 1 building 10 floor 2 appVII room 3 Moscow Russia 111524

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Test tubes are potentially hazardous medical waste to the spread of infectious diseases transmitted through blood and are class B medical waste — epidemiologically hazardous waste. Collection and disposal of used tubes should be carried out in accordance with the requirements of sanitary regulations of SanPin 2.1.7.2790-10 dated 12.12.2010 «Sanitation and epidemiological requirements for medical waste circulation» and MU 287-113 "Guidelines for disinfection, sterilizing cleaning and sterilization of medical products".

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REAL CARE LLC

Perovskayast. 1 building 10 floor 2 appVII room 3 Moscow Russia 111524

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
INSTRUCTION FOR USE

proved
Name: BAI QIANG
Position: SALES MANAGER

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Table 1

Name of the tube	Cap (plug) color
Vacuum tube VacPlus without filler	red

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Product is intended for use in collecting, storing and / or transporting venous blood samples for tests and / or other research, obtaining blood components.

Product is intended for use in the clinical diagnostic laboratories of medical institutions, blood transfusion stations and research institutions.

Purpose: Test tubes are used for research of serum. Test tubes are used for clinical, chemical and serological studies.

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INSTRUCTION FOR USE

- 1) Identify the patient, introduce himself, and explain the course and purpose of the procedure. Ensure that the patient has informed consent for the upcoming blood collection procedure. In the absence of such, clarify further actions with a doctor.
- 2) Offer the patient or help him to take a comfortable position: sitting or lying down.
- 3) Mark the tubes, indicating the patient's full name, department (in order to eliminate errors in the identification of the biomaterial sample). Using the digital code make appropriate notes in the logbook. (Be sure to put on the label on the tube, the name of the medical organization, the name and initials of the patient, his gender, age, date and time of blood collection and send the tube to the laboratory. Labels should always be filled in the presence of the patient.)
- 4) Treat your hands hygienically, dry them.
- 5) Treat your hands with an antiseptic. Do not dry, wait until the antiseptic is completely dry.
- 6) Wear medical gloves.
- 7) Select tubes corresponding to the declared tests or laboratory tests required by the patient, prepare a needle, holder, alcohol wipes, and a plaster.

Apply a tourniquet on a shirt or diaper 7-10 cm above the venipuncture site. Harness need to be no more than one minute. A longer vein squeeze time may affect the test results due to changes in the concentration of certain parameters in the blood.

Ask the patient to make a fist. It is impossible to set for a hand an exercise stress (vigorous squeezing and unclenching a fist"), since this may lead to changes in blood concentrations of some indicators. Choose venipuncture site. The most common are the ulnar and saphenous veins, but smaller and fuller veins of the dorsal surface of the wrist.

Performing the procedure:

- 1) Take the needle and remove the protective cap from it. If a double-sided needle is used, remove the protective cap in gray or white.
 - 2) Insert the needle into the needle holder and screw in until it stops.
 - 3) Disinfect the venipuncture site with a gauze napkin or swab moistened with an antiseptic solution, in a circular motion from the center to the periphery.
 - 4) Wait until the antiseptic solution is completely dry (30-60 s). You cannot wipe and blow the puncture site, so as not to bring microorganisms. It is also impossible to palpate the vein after disinfection. If difficulties arise during venipuncture and the vein is palpated again, this area must be disinfected again.
 - 5) Remove the cap on the other side of the needle.
 - 6) To clasp the patient's forearm with the left hand so that the thumb is 3-5 cm below the venipuncture site.
 - 7) Position the needle along the line with the vein bevel up and puncture the vein at an angle 15-30 ° to the skin.
 - 8) Insert the previously prepared tube into the needle holder until it stops and hold it until the blood stops flowing into the tube. The harness must be removed immediately after the start of blood flow into the tube. Make sure the patient opens his fist. The blood passes into the tube until it fully compensates for the vacuum created in it. If the blood does not flow, it means that the needle went through the vein - in this case, you need to slightly pull out the needle, but do not remove it until the blood goes to the tube. The accuracy of filling the tube is $\pm 10\%$ of the nominal volume.
 - 9) Remove the tube from the holder.
 - 10) Immediately after filling, the tube must be gently inverted to mix the sample with the filler: a tube without anticoagulants - 5-6 times; test tube with citrate - 3-4 times; a tube with heparin, EDTA and other additives - 8-10 times. The tubes cannot be shaken - it can cause foaming and hemolysis, as well as lead to mechanical lysis of red blood cells.
- If necessary, a number of other tubes are inserted into the needle holder to obtain the desired volume of blood for various studies. Re-insert the needle into the vein for this is not necessary.
- 11) After all the required tubes are filled, attach a dry sterile cloth to the venipuncture site and remove the needle.

III Termination of the procedure:

- 1) Place the used needle together with the disposable holder in a sharps container. In order to prevent contact with blood, it is forbidden to disassemble the needle and the holder in your hands!
- 2) Ensure that the patient does not have external bleeding in the venipuncture area.
- 3) Apply a pressure bandage on the arm or a bactericidal patch.
- 4) Handle gloved hands with a disinfectant.
- 5) Disinfect all used material. Remove gloves, place in a disinfection tank or a waterproof bag / container for disposal of class B waste.
- 6) Treat your hands hygienically, dry them.
- 7) Refine patient well-being.
- 8) Make an appropriate record of the results of the service in the medical records or arrange a referral.

...ge the delivery of test tubes with the obtained laboratory material to the laboratory
...e to put on the label on the tube, the name of the medical organization, the patient's name and
...s, his gender, age, date and time of blood collection and send the tube to the laboratory. Labels
... always be filled in the presence of the patient.

Transport labeled tubes in an upright position to appropriate laboratories in special containers
... with lids to be disinfected.

When taking blood from a peripheral vein for studies using various reagents, the following
sequence must be observed: first, blood is taken without anticoagulants, then in test tubes with
anticoagulants to avoid contamination.


Centrifuging


If necessary, the centrifugation procedure, make sure that the tubes are inserted into the rotor
in such a way that the cap does not rest on the walls of the centrifuge cup, otherwise the cap may
come off the tubes - a more stable gel barrier is formed in horizontal centrifuges than in centrifuges
with a fixed angle of inclination - high temperatures can have a negative effect on the physical
properties of the gel. The most complete separation of serum and plasma occurs at a temperature of
20°C - 22°C.

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Additive	Maximum time after blood sampling	How many times to flip	Recommended relative centrifugal acceleration	Recommended centrifugation time, min
No additive	24 hours	-	1500 g	10

SYMBOLS WITH DESCRIPTION

	Batch code
	Used by date
	Manufacturer
	Date of manufacturer
	Do not reuse
	Do not use if pack is damaged

	Sterilized by radiation
	Keep away from sunlight
	Keep dry
	Temperature limit 4-25°C
	Fragile
	Up
	Layer constraint
	CE marking
	In vitro diagnostic medical device

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Stopper-Butyl rubber

STERILIZATION

devices are sterile. Sterilization method – radiation.

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Questions about the quality of a medical product, manufactured by HEBEI XINLE SCI&TECH
D, China, please contact:
CARE LLC
ovskayast. 1 building 10 floor 2 appVII room 3 Moscow Russia 111524
Telephone: +7(495)234-57-31
E-mail: lab@int-med.ru

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CARE LLC

skayast. 1 building 10 floor 2 appVII room 3 Moscow Russia 111524

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of HEBEI XINLE SCI & TECH CO.,LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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Vacuum tube VacPlus with Clot-activator (SiO₂)
Vacuum tube VacPlus with Clot-activator (SiO₂) and separation gel
Vacuum tube VacPlus with Clot-activator (SiO₂) and G-Plus mid-located gel

RC № XXX от XX.XX.XXXX

TITLE. DESCRIPTION

1. Vacuum tube VacPlus with Clot-activator (SiO₂) (red). Blood volume: 1 ml, 2 ml, 3 ml, 4 ml, 5 ml, 6 ml, 7 ml, 8 ml, 9 ml, 10 ml. Dimensions: 13 mm x 75 mm, 13 mm x 100 mm, 16 mm x 100 mm.
2. Vacuum tube VacPlus with Clot-activator (SiO₂) and separation gel (yellow). Blood volume: 3 ml, 3.5 ml, 4 ml, 5 ml, 6 ml, 7 ml, 8 ml, 8.5 ml. Dimensions: 13 mm x 75 mm, 13 mm x 100 mm, 16 mm x 100 mm.
3. Vacuum tube VacPlus with Clot-activator (SiO₂) and G-Plus mid-located gel – (pink). Blood volume: 3 ml, 3.5 ml, 4 ml, 5 ml, 6 ml, 7 ml, 8 ml, 8.5 ml. Dimensions: 13 mm x 75 mm, 13 mm x 100 mm, 16 mm x 100 mm.

The product is intended for use in collecting, storing and / or transporting venous blood samples for analysis and / or other research, obtaining blood components.
The product is intended for use in the clinical diagnostic laboratories of medical institutions, blood transfusion stations and research institutions.

Type of sample - Venous blood sample.

Blood collection using a vacuum tube is possible with a needle inserted into a vein. Thanks to the vacuum created during manufacturing, a strictly dosed amount of venous blood enters the tube.

This product consists of a plastic tube, a rubber stopper, a protective cap, a label and additive: Clot activator SiO₂, Clot-activator SiO₂ with separation gel, Clot-activator (SiO₂) and G-Plus mid-located gel. Gel for separating serum from fibrin and blood cells, forming a strong barrier, ensuring the stability of the composition of the serum.

The tube is made of PET, which practically does not beat and is a clinically tested, safer alternative to glass. The stopper tightly closes the tube, regardless of the presence of a vacuum in it. The protective cap is designed as a plug cap for easy

connection and detachment, which eliminates the possibility of direct contact with the health worker's hands and blood samples that may cause sample contamination.

All tubes are sterile, equipped with a convenient label with a field for recording and double tear-off digital code. Each tube has a mark to accurately determine the level of filling (filling line) located on the label or on the tube itself, depending on the intended volume. Compliance with the required level of filling allows to achieve the correct ratio of blood and filler (additive). Test tubes are available in a wide range of volumes. Test tube sizes are suitable for most centrifuge models. The inner walls of the tubes are covered with fillers (additives) applied by micro dispersed spraying. The caps of the test tubes are color-coded in accordance with the international standard ISO 6710.

Table 1

Name of the tube	Cap (plug) color
Vacuum tube VacPlus with Clot-activator (SiO ₂)	Red
Vacuum tube VacPlus with Clot-activator (SiO ₂) and separation gel	Yellow
Vacuum tube VacPlus with Clot-activator (SiO ₂) and G-Plus mid-located gel	pink

Test tubes with clot activator. The inner walls of the tubes are coated with silica microparticles, which activate the process of clot formation when the tube is carefully turned over. The recommended clotting time is 10 to 30 minutes.

Tubes with clot activator and separation gel contain a separation gel located at the bottom of the tube. The density of this material is lower than that of a blood clot, but higher than that of serum. During centrifugation, the gel moves to the boundary between the clot and serum, where it forms a stable barrier. Thanks to this, serum can be taken directly from the tube. In gel tubes, serum remains stable for up to 48 hours under recommended storage conditions.

Attention! Gel tubes cannot be transported without prior centrifugation! Sudden movements can destroy the gel and cause its particles to enter the sample, as a result of which the analysis results will be erroneous. If transportation is required, then after centrifugation, it is necessary to keep the gel tubes in a rack at room temperature for at least 1 hour to form a strong gel barrier.

DESIGNATION AND INDICATIONS (purpose)

The product is intended for use in collecting, storing and / or transporting venous blood samples for analysis and / or other research, obtaining blood components.

The product is intended for use in the clinical diagnostic laboratories of medical institutions, blood transfusion stations and research institutions.

Purpose: Test tubes are used for research of serum. Test tubes are used for biochemical, immunological studies.

CONTRAINDICATIONS AND SIDE EFFECTS

There are no contraindications and side effects at this time.

PRECAUTIONS. WARNINGS. CAUTIONS

Precautions

- Handle all biological samples and sharp objects in accordance with the rules adopted at your medical facility
1. Discard used sharps into a special container for disposal.
 3. All liquid reagents are clear. Do not use tubes if the additives have changed color or there is a foreign substance in them.
 4. Do not use tubes after the expiration date.

Warnings:

1. Please check the directions for use and designation.
2. Check the tube for damage, contamination, leakage.
3. Verify that the tube specification and test parameter match.

Cautions and storage requirements (indicated on the label):

- Store at temperature from +4 ° C to 25 ° C
- Do not use if packaging is opened or damaged.
- Do not use after expiration date.

INSTRUCTION FOR USE

- 1) Identify the patient, introduce himself, and explain the course and purpose of the procedure. Ensure that the patient has informed consent for the upcoming blood collection procedure. In the absence of such, clarify further actions with a doctor.
- 2) Offer the patient or help him to take a comfortable position: sitting or lying down.
- 3) Mark the tubes, indicating the patient's full name, department (in order to eliminate errors in the identification of the biomaterial sample). Using the digital code make appropriate notes in the logbook.
(Be sure to put on the label on the tube, the name of the medical organization, the name and initials of the patient, his gender, age, date and time of blood collection and send the tube to the laboratory. Labels should always be filled in the presence of the patient.)
- 4) Treat your hands hygienically, dry them.
- 5) Treat your hands with an antiseptic. Do not dry, wait until the antiseptic is completely dry.
- 6) Wear medical gloves.
- 7) Select tubes corresponding to the declared tests or laboratory tests required by the patient, prepare a needle, holder, alcohol wipes, and a plaster.
- 8) Place a tourniquet on a shirt or diaper 7-10 cm above the venipuncture site. Harness need to impose no more than one minute. A longer vein squeeze time may affect the test results due to changes in the concentration of certain parameters in the blood.
- 9) Ask the patient to make a fist. It is impossible to set for a hand an exercise stress (vigorous "squeezing and unclenching a fist"), since this may lead to changes in blood concentrations of some indicators. Choose venipuncture site. The most common are the ulnar and saphenous veins, but smaller and fuller veins of the dorsal surface of the wrist.

Performing the procedure:

- 1) Take the needle and remove the protective cap from it. If a double-sided needle is used, remove the protective cap in gray or white.
- 2) Insert the needle into the needle holder and screw in until it stops.
- 3) Disinfect the venipuncture site with a gauze napkin or swab moistened with an antiseptic solution, in a circular motion from the center to the periphery.
- 4) Wait until the antiseptic solution is completely dry (30-60 s). You cannot wipe and blow the puncture site, so as not to bring microorganisms. It is also impossible to palpate the vein after disinfection. If difficulties arise during venipuncture and the vein is palpated again, this area must be disinfected again.

Remove the cap on the other side of the needle.

To clasp the patient's forearm with the left hand so that the thumb is 3-5 cm below the venipuncture site.

7) Position the needle along the line with the vein bevel up and puncture the vein at an angle 15-30 ° to the skin.

8) Insert the previously prepared tube into the needle holder until it stops and hold it until the blood stops flowing into the tube. The harness must be removed immediately after the start of blood flow into the tube. Make sure the patient opens his fist. The blood passes into the tube until it fully compensates for the vacuum created in it. If the blood does not flow, it means that the needle went through the vein - in this case, you need to slightly pull out the needle, but do not remove it until the blood goes to the tube. The accuracy of filling the tube is $\pm 10\%$ of the nominal volume.

9) Remove the tube from the holder.

10) Immediately after filling, the tube must be gently inverted to mix the sample with the filler: a tube without anticoagulants - 5-6 times; test tube with citrate - 3-4 times; a tube with heparin, EDTA and other additives - 8-10 times. The tubes cannot be shaken - it can cause foaming and hemolysis, as well as lead to mechanical lysis of red blood cells. If necessary, a number of other tubes are inserted into the needle holder to obtain the desired volume of blood for various studies. Re-insert the needle into the vein for this is not necessary.

11) After all the required tubes are filled, attach a dry sterile cloth to the venipuncture site and remove the needle.

III Termination of the procedure:

1) Place the used needle together with the disposable holder in a sharps container. In order to prevent contact with blood, it is forbidden to disassemble the needle and the holder in your hands!

2) Ensure that the patient does not have external bleeding in the venipuncture area.

3) Apply a pressure bandage on the arm or a bactericidal patch.

4) Handle gloved hands with a disinfectant.

5) Disinfect all used material. Remove gloves, place in a disinfection tank or a waterproof bag / container for disposal of class B waste.

6) Treat your hands hygienically, dry them.

7) Refine patient well-being.

8) Make an appropriate record of the results of the service in the medical records or arrange a referral.

9) Arrange the delivery of test tubes with the obtained laboratory material to the laboratory

Be sure to put on the label on the tube, the name of the medical organization, the patient's name and initials, his gender, age, date and time of blood collection and send the tube to the laboratory. Labels must always be filled in the presence of the patient.

Transport labeled tubes in an upright position to appropriate laboratories in special containers with lids to be disinfected.

When taking blood from a peripheral vein for studies using various reagents, the following sequence must be observed: first, blood is taken without anticoagulants, then in test tubes with anticoagulants to avoid contamination.

Centrifuging

If necessary, the centrifugation procedure, make sure that the tubes are inserted into the rotor in such a way that the cap does not rest on the walls of the centrifuge cup, otherwise the cap may come off the tubes - a more stable gel barrier is formed in


horizontal centrifuges than in centrifuges with a fixed angle of inclination - high temperatures can have a negative effect on the physical properties of the gel. The most complete separation of serum and plasma occurs at a temperature of 20°C - 22°C.


Test tubes with gel must be centrifuged no later than 2 hours after blood collection. When the gel barrier is already formed, the tubes should not be re-centrifuged. Re-centrifuging the tubes can result the destruction of the gel barrier and the entry of gel particles into the serum / plasma.

Preparation of samples for analysis

Additive	Maximum time after blood sampling	How many times to flip	Recommended relative centrifugal acceleration	Recommended centrifugation time, min
Clot activator	3 hours	5-6 times	1500 g	10
Clot activator + gel	3 hours before centrifugation and 48 hours after centrifugation	5-6 times	1800 g	10

SYMBOLS WITH DESCRIPTION

	Batch code
	Used by date
	Manufacturer
	Date of manufacturer
	Do not reuse

	Do not use if pack is damaged
STERILE R	Sterilized by radiation
	Temperature limit 4-25°C
	Keep dry
	Keep away from sunlight
	Fragile
	Up
	Layer constraint
CE	CE marking
IVD	In vitro diagnostic medical device

GUARANTEED BY MANUFACTURER VALUES OF MAIN PARAMETERS, CHARACTERISTICS (FEATURES) OF MEDICAL DEVICE.
 Values of the parameters, characteristics (features) of the medical device are presented in the table below. Where applicable, the permissible deviations from the nominal characteristics are $\pm 10\%$, unless otherwise agreed.

Name	Volume	Dimensions	Additive	Cap (plug) color	Quantity
Vacuum tube VacPlus with Clot-activator (SiO ₂) (red). Blood volume: 1 ml, 2 ml, 3 ml, 4 ml, 5 ml, 6 ml, 7 ml, 8 ml, 9 ml, 10 ml. Dimensions: 13 mm x 75 mm, 13 mm x 100 mm, 16 mm x 100 mm.	1 ml, 2 ml, 3 ml, 4 ml, 5 ml, 6 ml, 7 ml, 8 ml, 9 ml, 10 ml	13x75 mm, 13x100 mm, 16x100 mm	Clot-activator (SiO ₂)	Red	50, 100 pcs/rack
Vacuum tube VacPlus with Clot-activator (SiO ₂) and separation gel (yellow). Blood volume: 3 ml, 3.5 ml, 4 ml, 5 ml, 6 ml, 7 ml, 8 ml, 8.5 ml. Dimensions: 13 mm x 75 mm, 13 mm x 100 mm, 16 mm x 100 mm.	3 ml, 3.5 ml, 4 ml, 5 ml, 6 ml, 7 ml, 8 ml, 8.5 ml	13x75 mm, 13x100 mm, 16x100 mm	Clot-activator (SiO ₂) and separation gel (composition: Vinyl silicone oil 65%(CAS 68083-19-2), Gas phase silica 25%(CAS 14808-60-7), Dimethyl silica Gel 10%(CAS 112926-00-8))	Yellow	50, 100 pcs/rack
Vacuum tube VacPlus with Clot-activator (SiO ₂) and G-Plus mid-located gel – (pink). Blood volume: 3 ml, 3.5 ml, 4 ml, 5 ml, 6 ml, 7 ml, 8 ml, 8.5 ml. Dimensions: 13 mm x 75 mm, 13 mm x 100 mm, 16 mm x 100 mm.	3 ml, 3.5 ml, 4 ml, 5 ml, 6 ml, 7 ml, 8 ml, 8.5 ml	13x75 mm, 13x100 mm, 16x100 mm	Clot Activator (SiO ₂) and mid-located gel (composition: Vinyl silicone oil 65%(CAS 68083-19-2), Gas phase silica 25%(CAS 14808-60-7), Dimethyl silica Gel 10%(CAS 112926-00-8),)	Pink	50, 100 pcs/rack

MAINRAW MATERIALS

Tube - PET

Cap-PE

Stopper-Butyl rubber

STERILIZATION

The devices are sterile. Sterilization method – radiation.

PACKAGE

Test tubes are packed in racks of 100 pieces. or 50 pcs. and wrapped in shrink film (group packaging).

Test tubes are delivered together with the application instruction (1 instruction for group packaging).

OPERATION, STORAGE, TRANSPORTATION CONDITIONS. SHELF LIFE

Transportation of products without blood should be carried out by all types of covered transport in compliance with the conditions and requirements established for this type of transport, at the temperature from 4 to +25 ° C. During transportation to provide protection from the sun, rain and snow, do not mix with harmful or toxic substances. No special requirement stances.

Transportation

Temperature	+4 - +25°C
Atmospheric pressure	700 - 1060 hPa
Relative humidity	Not more 75%

Storage: Keep vacuum tubes for taking blood without blood at a temperature of 4-25 ° C in a dry ventilated room. Keep away from harmful or poisonous substances.

Attention! Avoid direct sunlight. Exceeding the recommended storage temperature may lead to a loss of quality of the mini tube (for example, evaporation or staining of fillers (additives)).

Storage

Temperature	+4 - +25°C
Atmospheric pressure	700 - 1060 hPa
Relative humidity	Not more 75%

Operation

Temperature	15-30° C
Atmospheric pressure	700 - 1060 hPa
Relative humidity	35% - 75%

Shelf life

Shelf life is 24 months

WARRANTY OBLIGATION

All questions related to the operation and quality of the product can be addressed to the Authorized representative of the manufacturer in the territory of the Russian Federation. The manufacturer guarantees that the product meets the requirements of the standards and documentation of the manufacturer under the conditions of transportation, storage, operation.

The manufacturer shall not be liable for any damage or damage caused by the use of the product until proven guilty.

The manufacturer shall not be liable for any damage or damage caused by the use of the product after the expiry date indicated on the packaging.

For questions about the quality of a medical product, manufactured by HEBEI XINLE SCI&TECH CO.,LTD, China, please contact:

REAL CARE LLC

Perovskaya st. 1 building 10 floor 2 appVII room 3 Moscow Russia 111524

Telephone: +7(495)234-57-31

E-mail:lab@int-med.ru

POTENTIAL CONSUMER

A specialist with a diploma in the established form on graduating from a secondary professional medical educational institution in the specialties: "medical care", "nursing", "obstetric care".

There are skills to perform this simple medical service.

Patient without age restrictions.

Execution conditions

Outpatient.

Stationary

THE MAINTENANCE AND REPAIR OF MEDICAL DEVICES

Not applicable.

DISPOSAL OR DESTRUCTION OF MEDICAL DEVICES.

Used test tubes are potentially hazardous medical waste to the spread of infectious diseases transmitted through blood and are class B medical waste — epidemiologically hazardous waste. Disinfection and disposal of used tubes should be carried out in accordance with the requirements of sanitary regulations of SanPin 2.1.7.2790-10 dated 12.12.2010 «Sanitation and epidemiological requirements for medical waste circulation» and MU 287-113 "Guidelines for disinfection, presterilizing cleaning and sterilization of medical products".

Unused tubes, including expired ones, should be disposed of as Class A waste and can be disposed of as household waste.

Adhere to national regulations when disposing or recycling the product, its components and its packaging.

5

IFUvac(SiO₂, SiO₂+gel, SiO₂+G-Plus gel): 01 - 2020

A blood sample is a hidden source of contamination. Therefore, all components of the equipment, consumables and waste substances that come into contact with the sample should be used in accordance with biological safety rules.

Manufacturer:

Hebei Xinle Sci&Tech Co., Ltd China

No. 189, Nanhuan Road 050700 Xinle City, Hebei Province People's Republic of China

+86 311 8858 1771

hbxinle001@hbxinle.net

Manufacturing site:

Hebei Xinle Sci&Tech Co., Ltd, China

No. 189, Nanhuan Road 050700 Xinle City, Hebei Province People's Republic of China

AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION:

REAL CARE LLC

Perovskaya st. 1 building 10 floor 2 appVII room 3 Moscow Russia 111524

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证明书查询 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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Vacuum tube VacPlus with EDTA K3
Vacuum tube VacPlus with EDTA K3 and separation gel
Vacuum tube VacPlus with EDTA K2
Vacuum tube VacPlus with EDTA K2 and separation gel

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This product consists of a plastic tube, a rubber stopper, a protective cap, a label and additive: EDTA K2, EDTA K3, EDTA K2 and gel, EDTA K3 and gel.

The tube is made of PET, which practically does not break and is a clinically tested, safer alternative to glass. The stopper tightly closes the tube, regardless of the presence of a vacuum in it. The protective cap is designed as a plug cap for easy connection and detachment, which eliminates the possibility of direct contact with the health worker's hands and blood samples that may cause sample contamination.

All tubes are sterile, equipped with a convenient label with a field for recording with double tear-off digital code. Each tube has a mark to accurately determine the level of filling (filling line) located on the label or on the tube itself, depending on the intended volume. Compliance with the required level of filling allows to achieve the correct ratio of blood and filler (additive). Test tubes are available in a

range of volumes. Test tube sizes are suitable for most centrifuge models. The inner walls of the tubes are covered with fillers (additives) applied by micro dispersed spraying. The caps of the test tubes are color-coded in accordance with the international standard ISO 6710.

Table 1

Designation of the tube	Cap (plug) color
Vacuum tube VacPlus with EDTA K3	lilac
Vacuum tube VacPlus with EDTA K3 and separation gel	
Vacuum tube VacPlus with EDTA K2	
Vacuum tube VacPlus with EDTA K2 and separation gel	

The inner walls of the tubes are coated with microdrops of the anticoagulant EDTA K2 (2-substituted potassium salt of ethylenediaminetetraacetic acid) or EDTA K3 (3-substituted potassium salt of ethylenediaminetetraacetic acid). EDTA prevents blood coagulation by binding calcium ions and does not affect hematological parameters. A blood smear must be done within 3 hours after sampling, the remaining studies are carried out for 6 hours at room temperature or for 24 hours at a temperature of 0-2 ° C. The red blood cells, white blood cells and platelets in the blood sample with the EDTA anticoagulant are stable up to 24 hours.

Attention! If transportation of samples is required, tubes with EDTA K2 should be used, since EDTA K3 exhibits slightly lower properties as an anticoagulant and much more affects the morphology of cells during long-term storage (more than 4 hours).

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Purpose: Test tubes are used for hematological, immunochemical, virological and qualitative and quantitative studies of whole blood.

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Warnings:

1. Please check the directions for use and designation.
2. Check the tube for damage, contamination, leakage.

Verify that the tube specification and test parameter match.

Conditions and storage requirements (indicated on the label):

Store at temperature from +4 ° C to 25 ° C

Do not use if packaging is opened or damaged.

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Performing the procedure:

- 1) Take the needle and remove the protective cap from it. If a double-sided needle is used, remove the protective cap in gray or white.
- 2) Insert the needle into the needle holder and screw in until it stops.
- 3) Disinfect the venipuncture site with a gauze napkin or swab moistened with an antiseptic solution, in a circular motion from the center to the periphery.
- 4) Wait until the antiseptic solution is completely dry (30-60 s). You cannot wipe and blow the puncture site, so as not to bring microorganisms. It is also impossible to palpate the vein after disinfection. If difficulties arise during venipuncture and the vein is palpated again, this area must be disinfected again.
- 5) Remove the cap on the other side of the needle.
- 6) To clasp the patient's forearm with the left hand so that the thumb is 3-5 cm below the venipuncture site.
- 7) Position the needle along the line with the vein bevel up and puncture the vein at an angle 15-30 ° to the skin.
- 8) Insert the previously prepared tube into the needle holder until it stops and hold it until the blood stops flowing into the tube. The harness must be removed immediately after the start of blood flow into the tube. Make sure the patient opens his fist. The blood passes into the tube until it fully compensates for the vacuum created in it. If the blood does not flow, it means that the needle went

in the vein - in this case, you need to slightly pull out the needle, but do not remove it until the blood goes to the tube. The accuracy of filling the tube is $\pm 10\%$ of the nominal volume.

Remove the tube from the holder.

Immediately after filling, the tube must be gently inverted to mix the sample with the filler: a tube without anticoagulants - 5-6 times; test tube with citrate - 3-4 times; a tube with heparin, EDTA and other additives - 8-10 times. The tubes cannot be shaken - it can cause foaming and hemolysis, as well as lead to mechanical lysis of red blood cells.

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III Termination of the procedure:

- 1) Place the used needle together with the disposable holder in a sharps container. In order to prevent contact with blood, it is forbidden to disassemble the needle and the holder in your hands!
 - 2) Ensure that the patient does not have external bleeding in the venipuncture area.
 - 3) Apply a pressure bandage on the arm or a bactericidal patch.
 - 4) Handle gloved hands with a disinfectant.
 - 5) Disinfect all used material. Remove gloves, place in a disinfection tank or a waterproof bag / container for disposal of class B waste.
 - 6) Treat your hands hygienically, dry them.
 - 7) Refine patient well-being.
 - 8) Make an appropriate record of the results of the service in the medical records or arrange a referral.
 - 9) Arrange the delivery of test tubes with the obtained laboratory material to the laboratory
- Be sure to put on the label on the tube, the name of the medical organization, the patient's name and initials, his gender, age, date and time of blood collection and send the tube to the laboratory. Labels must always be filled in the presence of the patient.

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Additive	Maximum time after blood sampling	How many times	Recommended relative centrifugal	Recommended centrifugation time,
----------	-----------------------------------	----------------	----------------------------------	----------------------------------

		to flip	acceleration	min
EDTA K2	24 hours Test tubes with gel must be centrifuged no later than 2 hours after blood collection	8-10 times	1300 g	10
EDTA K3	24 hours Test tubes with gel must be centrifuged no later than 2 hours after blood collection	8-10 times	1300 g	10

SYMBOLS WITH DESCRIPTION

	Batch code
	Used by date
	Manufacturer
	Date of manufacturer
	Do not reuse
	Do not use if pack is damaged
	Sterilized by radiation
	Keep away from sunlight
	Keep dry
	Temperature limit 4-25°C
	Fragile

	Up
	Layer constraint
	CE marking
	In vitro diagnostic medical device

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