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
The screenshot shows the 'Ввести рецепт' (Enter recipe) window. It contains several input fields and buttons. The 'Напряж.' (Voltage) is set to 50 [кВ], 'Ток' (Current) to 40 [μA], and 'Тип' (Type) to 'Сферич' (Spherical). The 'Applicator' is set to '3.0' and 'ID пациен.' (Patient ID) to 'ID1234'. The 'Лечение №' (Treatment No.) is '411'. The 'Предпис. доза' (Prescribed dose) is '20.00 [Gy]' and 'Время лечения' (Treatment time) is '0:23:51 [ч:м:с]'. The 'Глуб. возд-вия' (Depth of action) is '0.00 мм от поверхности' and 'Мощ. дозы при глуб.' (Power of dose at depth) is '0.838 [Gy/мин]'. The 'Рецепт готов' (Recipe ready) button is visible.

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- Armoogum KS, Parry JM, Souliman SK, Sutton DG, Mackay CD. Functional intercomparison of intraoperative radiotherapy equipment – Photon Radiosurgery System. *Radiat Oncol* 2007;2:11.
- Armoogum K, Parry J, Mackay C, Thompson A, Vaidya J, Ackland C, Gardner J, Munro A, Eljamel S: Implementation and experiences of an intraoperative radiotherapy service. *Clin Oncol (R Coll Radiol)* 2007;19(3 Suppl):S8-9.
- Ebert MA, Carruthers B: Dosimetric Characteristics of a Low-kV intra-operative X-ray Source: Implications for use in a Clinical Trial for Treatment of Low-risk Breast Cancer. *Med Phys* 2003;30:2424-31.
- Ebert MA, Carruthers B, Lanzon PJ, Haworth A, Clarke J, Caswell NM, Siddiqui SA: Dosimetry of a Low-kV Intra-operative X-ray Source using Basic Analytical Beam Models. *Australas Phys Eng Sci Med* 2002;25:119-23.
- Clarke J, Lanzon P, Haworth A, Ebert M, Caswell N, Tuchnya T: The PeC Intrabeam Miniature X-Ray Device – Physics Implementation. *Proceedings Intercontinental Workshop and Symposium on IORT using Mobile Systems, Sydney, 2002.*
- Biggs PJ, Beatty J, Yasuda T: Distance-dose Curve for a Miniature X-ray Tube for Stereotactic Radiosurgery using an Optimized Aperture with a Parallel-Plate Ionization Chamber. *Med Phys* 1999;26:2550-54.
- Yasuda T, Beatty J, Biggs PJ, Gall K: Two-dimensional Dose Distribution of a Miniature X-ray Device for Stereotactic Radiosurgery. *Med Phys* 1998;25:1212 – 1216.

- Hakim R, Zervas NT, Hakim F, Butler WE, Beatty J, Yanch JC, Biggs PJ, Gall KP, Sliski AP: Initial Characterization of the Dosimetry and Radiobiology of a Device for Administering Interstitial Stereotactic Radiosurgery. *Neurosurgery* 1997;40:510-7.
- Yanch JC, Harte KJ: Monte Carlo Simulation of a Miniature, Radiosurgery X-ray tube using the IST 3.0 Coupled Electron-Photon Transport Code. *Med Phys* 1996;23:1551-8.
- Dinsmore M, Harte KJ, Sliski AP, Smith DO, Nomikos PM, Dalterio MJ, Boom AJ, Leonard WF, Oettinger PE, Yanch JC: A new Miniature X-ray Source for Interstitial Radiosurgery: Device Description. *Med Phys* 1996;23:45-52.
- Beatty J, Biggs PJ, Gall K, Okunieff P, Pardo FS, Harte KJ, Dalterio MJ, Sliski AP: A new Miniature X-ray Device for Interstitial Radiosurgery: Dosimetry. *Med Phys* 1996;23:53-62.
- Biggs DS, Thomson ES: Radiation Properties of a miniature X-ray device for radiosurgery. *Br J Radiol* 1996;69:544-7

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- Herskind C, Kraus-Tiefenbacher U, Keller A, Fleckenstein K, Wenz F: Intraoperative Boost in mammary cancer using 50 kV X-Rays from INTRABEAM: Radiobiological modeling of late reaction. *Rev. Cancer* 2008; 22. No Extra:6.
- Belletti B et al., Targeted Intraoperative Radiotherapy Impairs the Stimulation of Breast Cancer Cell Proliferation and Invasion Caused by Surgical Wounding, *Clin Cancer Res* 1325 2008;14(5)
- Gardner J, Ackland CY, Thompson AM, Valdyia SJ, Dewar JA, Adamson DJ, Armoogum K, Mackay CD, Munro AJ: On the Practicalities of Translational Radiobiology: Lessons from IORT (Intra-operative Radiotherapy). *Clin Oncol (R Coll Radiol)* 2007;19(3 Suppl):S29.
- Herskind C, Schalla S, Hahn EW, Hver KH, Wenz F: Influence of different dose rates on cell recovery and RBE at different spatial positions during protracted conformal radiotherapy. *Radiat Prot Dosimetry* 2006;122:498-505.

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- Vaidya J et al. Targeted intraoperative radiotherapy (TARGIT) yields very low recurrence rates when given as a boost, *Int. J. Radiation Oncology Biol. Phys.*, Vol. 66, No. 5, pp. 1335–1338, 2006
- Vaidya J et al. Targeted intraoperative radiotherapy (TARGIT) Boost after breast conserving surgery results in a remarkably low recurrence rates in a standard risk population: 5 year results, *Proc. San Antonio Breast Cancer Conference* Dec 2009
- Wenz F et al. Intraoperative radiotherapy (IORT) as a boost during breast conserving surgery (BCS) using low kV X-rays - the first five years of experience with a novel approach, *Int. J. Radiation Oncology Biol. Phys.* (Article 18522)
- Kraus-Tiefenbacher U et al. Postoperative seroma formation after intraoperative radiotherapy using low-kilovoltage X-rays given during breast conserving surgery, *Int. J. Radiation Oncology Biol. Phys.*, 2009 Oct 14.
- Kraus-Tiefenbacher U et al. Long term toxicity of an intraoperative radiotherapy boost using low energy X-rays during breast conserving surgery, *Int. J. Radiation Oncology Biol. Phys.*, Vol. 66, No. 2, pp. 377-381, 2006
- Keshtgar M et al. Breast Cosmetic outcome 2 and 3 years after intraoperative radiotherapy compared with external beam radiotherapy for early breast cancer: An objective assessment of patients from a randomized controlled trial, *Breast Cancer symposium* 2009; A 294
- Kraus-Tiefenbacher U et al. Intraoperative Radiotherapy (IORT) as a Boost in Patients with Early-Stage Breast Cancer – Acute Toxicity, *Onkologie*. 2003 Dec;26(6):596-8
- Keshtgar M et al. Intraoperative Radiotherapy for Breast Cancer in patients where external beam radiation was not possible: Results after 3 years of follow up, *Proc. St Gallen Breast Cancer Conference* Mar 2009
- Kraus-Tiefenbacher U et al. Intraoperative radiotherapy (IORT) is an option for patients with localized breast recurrences after previous external-beam radiotherapy, *BMC Cancer* 2007, 7:178

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- Eljamel S et al. Intraoperative Radiotherapy using mobile miniature X-Ray Source in malignant brain tumors, Proc. American Association of Neurosurgery (AANS) 2008
- Kalapurakal JA, Goldmann S, Stelipflug W, Curran J, Sathaseelan V, Marymont MH, Tomita T: Phase I study of intraoperative radiotherapy with photon radiosurgery system in children with recurrent brain tumours: preliminary report of first dose level (10 Gy). Int J Rad Oncol Biol Phys 2006;65:800-8.
- Curry W et al. Stereotactic interstitial radiosurgery for cerebral metastases, J Neurosurg 103:630–635, 2005

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


ZEISS

INTRABEAM® System PRS 500 with XRS 4 Software-Version 2.2

Instructions for Use for the Core System

G-30-1601-en

Version 11.1

2016-09-01

ZEISS

- About this manual* The user manual is part of the delivery package.
- Carefully read it before using the system.
 - Keep it with the system.
 - Store it for the entire service life of the system.
 - Pass it on to subsequent owners or users of the system.
- Scope* These Instructions for Use are for the INTRABEAM® Core System PRS 500 with XRS 4 and Software Version 2.2, with Order Number 304534-0000-500.
- Trademark* INTRABEAM is a registered trademark of Carl Zeiss Meditec AG or of another company of the ZEISS Group in Germany and/or other countries.
- Registered trademark* UNIDOS® is a registered trademark of PTW.
- Manufacturer*
- | | |
|----------------------------|-----------------------------------|
| Carl Zeiss Meditec AG | |
| Goeschwitzer Strasse 51-52 | Fax: + 49 (0) 73 64 - 20 48 23 |
| 07745 Jena | Email: surgical@meditec.zeiss.com |
| Germany | Internet: www.meditec.zeiss.com |
- Original language* These Instructions for Use were originally drafted and published in German.
- Copyright* © Carl Zeiss Meditec AG 2016
All rights reserved. Reprints and reproductions, including excerpts, are not permissible. However, these Instructions for Use may be digitized and archived for internal use. Access to this data by third parties is not permissible and must be prevented.
- Change note* Subject to change in design and scope of delivery and as a result of ongoing technical development.

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Predecessor system	Previous INTRABEAM System PRS 400
XRS	X-ray source
PRS	Control console
ERM	External Radiation Monitor
ETA	ERM Test Adapter
EU	European Union
IC	Ionization Chamber
IRM	Internal radiation monitor
PAICH	Probe adjuster / ionization chamber holder
PDA	Photo Diode Array
SQA	System Quality assurance
EOT	End of Treatment
RTD	Resistance Temperature Device
EMC	Electromagnetic compatibility
PTB	Physikalisch Technische Bundesanstalt
NIST	U.S. National Institute of Standards and Technology
USA	United States of America
ABC-Error	Automatic Beam Centering Error)

Designations and abbreviations

Product name (English)	Designation in these Instructions for Use	
INTRABEAM system PRS 500 with XRS 4	INTRABEAM system	Currently: INTRABEAM System PRS 500 with XRS 4 and Software Version 2.2
INTRABEAM core system	INTRABEAM Core system	The elementary components of the INTRABEAM system form the core system: Control unit, user terminal, x-ray source and the components for checking: V-block, X-block, PDA, PAICH, ion- ization chamber and dosimeter.
Control unit PRS 500	Control console	Currently: PRS 500
XRS X-ray source		XRS = X-Ray Source Currently: XRS 4
Dosimeter	Dosimeter	Currently: Dosimeter UNIDOS E from PTW
PRS 500 Terminal	User terminal	Currently: PRS 500 user terminal
NC32 INTRABEAM Floor Stand	INTRABEAM stand	Currently: Floor Stand NC32 INTRABEAM
INTRABEAM drape	Drape	Steril drape for INTRABEAM stand
INTRABEAM cart	System cart	
INTRABEAM radia- tion shields	Shields	Foil for Shielding against X-radia- tion

Contents

	Abbreviations in these Instructions for Use	3
	Designations and abbreviations	4
Functions at a glance		
	INTRABEAM Core System	12
	Components for checking the system	14
	Accessories for checking the system	14
	Accessories of the INTRABEAM system	16
Safety measures		
	Key to symbols	21
	Hazard symbols	21
	Information symbols	21
	Target group	22
	Field of application	22
	Intended use	22
	Normal use	22
	Major performance characteristics	23
	Clinical experience	23
	Tips for owner/operators and users	24
	Use of the INTRABEAM Core System	30
	Data protection and information security	31
	Archiving of data	32
	Safety devices of the INTRABEAM Core System	33
	Calibrations	35
	Safety during tests and verification	35

Safety during treatment.....	36
Electrical features of the INTRABEAM Core System	37
Fire, anesthetics and oxygen	38
EMC (Electromagnetic compatibility)	39
Cleaning, disinfection, sterilization	39
Warning and information labels	43
PRS 500 Control Console rating plate.....	43
Warning label	45
Symbols and inscriptions laterally and at the bottom.....	46
Label on the XRS 4 X-ray source	47
Rating label of the PRS 500 user terminal.....	48

System overview

System description.....	50
Main components of the INTRABEAM Core System	50
Verification components	50
Components for the protection of the INTRABEAM Core System	51
Accessories of the INTRABEAM system	51
Principal components	53
PRS 500 control console	53
User terminal	58
XRS 4 X-ray radiation source and IRM internal X-ray radiation monitor	64
Protection components.....	66
X-block	66
V-block guide	67
Probe protector	68
Verification devices.....	70
Probe adjuster/ionization chamber holder PAICH	70
PDA photodiode array	72

Ionization chamber (IC).....	72
Cables - color coding.....	74
UNIDOS E dosimeter.....	75
Water phantom, optional verification accessories	76
Transport and storage container	78
Accessories for the INTRABEAM core system	80
INTRABEAM cart	80
NC32 INTRABEAM Floor Stand	81

Preparations for use

Requirements for system operation.....	85
Staff requirements.....	85
Space requirements.....	85
Electrical requirements.....	86
Visual system check.....	86
Positioning the INTRABEAM cart	87
Layout of the treatment facility.....	89
Connecting	92
Workflow	94
Powering up the INTRABEAM Core System.....	96
Working mode selection.....	98
Application software - main window	100
Menu bar	101
Tab bar.....	105
Data import and export	106
Equipment Inventory	112
"XRS X-Ray Sources" group	114
Recalibration of XRS X-ray sources.....	120

"QA Instruments" group	123
"Applicators" group	127
System Quality Assurance	130
Video instruction, video window	136
XRS probe adjustment using PAICH	138
Dynamic Offsets (adjusting the beam deflection)	142
PDA Source Check (isotropy)	146
PAICH Output Check (dose rate)	153
Only for predecessor systems: checking the ERM external radiation monitor. 161	
Entering treatment parameters	163
Checking the treatment parameters	173

Operation

Checklist for INTRABEAM Core System	178
Conducting treatment	180
Troubleshooting	191
Documenting the number of applicator sterilizations	196
Review treatment, displaying treatment parameters	199

Care and maintenance

Calibration regulations	209
Calibration	210
Date and time synchronization	212
Determining the deflection sensitivity	214
Determining the linearity of the IRM internal radiation monitor	218

Calibrating the air pressure sensor	220
Calibrating the temperature sensor	223
For predecessor systems only: Dilator Adjustment	226
Troubleshooting table	227
Changing the line fuse	233
Resterilization	234
Cleaning and disinfecting components	235
Cleaning the touchscreen	237
Cleaning the equipment trolley(INTRABEAM Cart)	238
Disinfecting user terminal and cart	238
Handling, transport and shipment	239
Disposal	240
Ordering data	241
Spare parts	242
Accessories	243
 System data	
Technical data of the XRS 4 x-ray source	246
Technical data for UNIDOS E	247
Technical data of the PRS 500 control console	248
Technical data of the user terminal	249
External interlock and X-radiation indicator	250
External X-ray interlock	250
External X-ray display	251
EMC (electromagnetic compatibility)	252
Electromagnetic radiation disturbance	254
Electromagnetic immunity	255

Recommended safety distances 258

Environmental conditions259

Regulatory information.....260

Regulatory information on the PRS 500 controller..... 260

Regulatory Information about the User Terminal..... 261

Supplementary documents

Applicable user manuals264

Equipment checklist for the operating room265

Literature.....266

Radiotherapy 266


Index

Functions at a glance

INTRABEAM Core System	12
Components for checking the system	14
Accessories for checking the system	14
Accessories of the INTRABEAM system	16

INTRABEAM Core System

PRS 500 control unit	Page 53
1 Emergency OFF switch on the control unit	Page 54
2 PAUSE button	Page 55
3 Speaker	Page 56
4 Power switch, line fuse and power inlet socket	Page 56
5 Potential equalization connector	Page 56
6 RS232 serial interface port	Page 92
/ PDA/PAICH connector (yellow)	Page 92
8 ERM connector (blue) - only for predecessor systems	Page 92
9 XRS connector (green)	Page 92
10 Connector for external interlock (gray)	Page 56
XRS 4 X-ray source	Page 64
11 XRS 4 x-ray source with integrated IRM x-ray monitor and with probe guard (11.1)	
Application software user interface	Page 96
Tabs for the activation of the following screen masks:	
12 Equipment Inventory	Page 112
13 System Quality Assurance	Page 130
14 "Enter prescription"	Page 163
15 "Verify treatment parameters"	Page 173
16 "Treatment"	Page 180
17 "Review Treatment"	Page 199
18 Calibration	Page 210
19 Heartbeat display	Page 100
PRS 500 user terminal	Page 58
20 Monitor (touchscreen) with computer	
21 Access to DVD drive and USB port	
22 Power ON key of user terminal	Page 96
23 Membrane keyboard with mouse pad	
24 Electrical connections on the underside of the device	


Components for checking the system

1	V-block guide (The illustration shows the V-block guide with the XRS X-ray source and PDA photo diode array.)	Page 67
2	X-block	Page 66
3	PDA photo diode array	Page 72
4	PAICH Special centering device and Ionization chamber holder	Page 70
5		Page 72
6	UNIDOS E dosimeter	Page 75

Accessories for checking the system

7	INTRABEAM Water Phantom	Page 76
---	-------------------------	---------



Accessories of the INTRABEAM system

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
Information on the accessories of the INTRABEAM system can be found in the corresponding separate Instructions for Use. Overview, see Page 26.

- | | |
|---|--|
| 1 | INTRABEAM Spherical Applicators (The illustration shows the applicator set in the sieve container) |
| 2 | INTRABEAM Flat Applicator Sets |
| 3 | INTRABEAM Surface Applicator Sets |
| 4 | INTRABEAM Needle Applicator |
| 5 | Radiation Shield, flat |
| 6 | NC32 INTRABEAM Floor Stand (INTRABEAM stand) |
| 7 | INTRABEAM Cart (equipment cart) |
| 8 | INTRABEAM Drapes |

(Blank page, for your notes . . .)

Safety measures

Key to symbols	21
Hazard symbols	21
Information symbols	21
Target group	22
Field of application	22
Intended use	22
Normal use	22
Major performance characteristics	23
Clinical experience	23
Tips for owner/operators and users	24
Use of the INTRABEAM Core System	30
Data protection and information security	31
Archiving of data	32
Safety devices of the INTRABEAM Core System	33
Calibrations	35
Safety during tests and verification	35
Safety during treatment	36
Electrical features of the INTRABEAM Core System	37
Fire, anesthetics and oxygen	38
EMC (Electromagnetic compatibility)	39
Cleaning, disinfection, sterilization	39
Warning and information labels	43
PRS 500 Control Console rating plate	43
Warning label	45
Symbols and inscriptions laterally and at the bottom	46

Label on the XRS 4 X-ray source	47
Rating label of the PRS 500 user terminal.....	48

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Additional information and tips

Target group

These Instructions for Use are intended for qualified specialist staff who have been trained and briefed in preparing and handling the medical device on the basis of these Instructions for Use.

Certain fields of application of the INTRABEAM system may require advanced qualification of users in the particular special field.

Field of application

Intended use

The "PRS 500 INTRABEAM System with XRS 4" is a system for radiation therapy treatment.

Normal use

INTRABEAM System PRS 500 with XRS 4 (INTRABEAM System) is a system for radiation therapy using x-ray beams.

The INTRABEAM system is a system for the minimal-invasive, highly focused treatment of selected lesions, on the one hand, and a system for intra-operative radiation therapy of the tumor bed within the entire body, on the other hand.

The core piece of the INTRABEAM System is the XRS miniature x-ray source XRS 4, equipped with an x-ray probe 10 cm in length and 3.2 mm in diameter. This x-ray probe delivers a high dose of low energy radiation.

In addition, electrons along the x-ray probe are accelerated using a maximum acceleration voltage of 50 kV. At the tip of the e-x-ray probe, they encounter a target object, and this process creates the x-ray beam. The effective emitted X-ray energy is approx. 20 keV or less. The dose quickly drops in soft tissue because of the low energy of the x-ray ($\sim 1/r^3$). This way, the treated lesion can shield the healthy tissue against injuries.

An electronic control unit acts as the voltage supply of the XRS X-ray source.

Various applicators are available, which can be connected to the probe tip of the XRS x-ray source in order to be put on a lesion for local radiation or introduced into the bed of a removed tumor.

The INTRABEAM System can be configured in such a way that it is possible to radiate defined lesions with interstitial, intra-operative and intra-cavity processes and to treat surfaces.

The NC32 INTRABEAM Floor Stand is available for the reliable positioning of the XRS x-ray source when an applicator is mounted. An INTRABEAM Drapes allows the XRS x-ray source and the stand arm to be covered and sterile.

The INTRABEAM Cart services you as an operating station for your INTRABEAM Core System and as a storage station for the components used to check the system and for accessories.

The INTRABEAM Water Phantom can be used to check the depth dose curve and the isotropy of your XRS x-ray source.

ZEISS will provide you with routs for quality assurance and for systematically recurring calibration work.

The INTRABEAM Core System is not intended for the archiving of data.

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The major performance characteristic of the INTRABEAM System is the delivery of defined, isotropic and constant x-ray radiation.

Clinical experience

Clinical trials of the INTRABEAM Core System began in December 1992. Results achieved in patients successfully demonstrate that the INTRABEAM core system is safe and effective for the irradiation of tumors.

Clinical experience to date includes use in open surgical resection of tumors in which the XRS X-ray source with ZEISS spherical applicators is used for intraoperative irradiation of the tumor bed following tumor resection.

Please refer to the bibliography listed in the chapter "Supplementary documents" for further information.

Tips for owner/operators and users

Responsibilities of the owner/ operator

The INTRABEAM System may only be used by qualified professionals, who have been trained how to operate the INTRABEAM System. The INTRABEAM system shall not be operated or used by other persons.

Every institution that operates an INTRABEAM System must appoint suitable and qualified expert staff who have been trained to operate the INTRABEAM System; this/those person(s) is/are responsible for the tasks involving the INTRABEAM System.

The manufacturer - Carl Zeiss Meditec AG - offers training courses on the INTRABEAM system.

The INTRABEAM system is restricted to purchase by or on the order of physicians or other authorized medical practitioners.

- Keep the Instructions for Use handy at all times for the operating staff.

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Precondition for additional equipment

Additional devices that are connected to electronic medical devices must comply with the relevant standards and guidelines. Furthermore, any configurations must comply with normative requirements for medical systems (see IEC 60601-1-1:2000 or

IEC 60601-1:2005, chapter 16). Anyone connecting additional equipment to medical electrical devices is a system configurer and as such responsible for compliance of the system with the normative requirements for systems.

Please note that local legislation takes priority over the above-mentioned normative requirements. Please contact ZEISS Service if you have any questions.

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Ambient conditions

- The INTRABEAM system should not be operated:
 - in areas under threat of explosion;
 - at a distance of less than 25 cm to any flammable anesthetics or volatile solvents such as alcohol, benzine or similar.
- Do not use or store components of the INTRABEAM system in damp spaces.

- Prevent dripping, surging or spraying water in the vicinity of components of the INTRABEAM system.
- Please ensure that no liquid is permitted to penetrate into the components of the INTRABEAM System.
- Ensure that the installation conditions and the use of the device meet surgical requirements:
 - Low vibration
 - Clean environment
 - Avoidance of extreme mechanical loading
- Observe the environmental conditions for the operation, transportation and storage in the "Device details" section, Page 259.

Electromagnetic compatibility (EMC)

- Please note the information on electromagnetic compatibility (EMC) in the chapter "System data", Page 252.

The electrically operated components of the INTRABEAM system comply with

- the EMC requirements of IEC 60601-1-2.
- Class B requirements for electromagnetic interference (practice area).
- It cannot be ruled out, however, that interferences may occur in HF receivers (e.g. TV or radio sets) in the vicinity. If such interference is noticed, please inform the ZEISS Service.
- When operating the system please observe the following precautions regarding EMC:
 - Only use accessory equipment approved by ZEISS for this system.
 - Do not use any portable or mobile RF communication equipment near the INTRABEAM system as it is not possible to exclude the possibility that the function of the device will be affected.

Warranty and Liability

Warranty and liability depend on the applicable contractual stipulations. The manufacturer will not be held liable for any damages, which are the result of unauthorized tampering with the components of the INTRABEAM System. Furthermore, this will forfeit any rights to claim under warranty.

- Do not undertake any modifications to the components of the INTRABEAM System without the manufacturer's permission.
- If a component is modified after this has been discussed with the manufacturer, further suitable inspections and tests will need to be performed to safeguard its safe continued use.

- Modifications and repairs of the components of the INTRABEAM System and the devices which are operated along with this INTRABEAM System may only be carried out by ZEISS Service or by authorized professionals.

Requirements from the user

The correct operation of the components of the INTRABEAM system is indispensable for the safe operation of the INTRABEAM System. Familiarize yourself thoroughly with the contents of the Instructions for Use prior to putting the INTRABEAM system into operation. Also observe the Instructions for Use of the other devices and the accessory components. Supplementary information is available from our Service department or from authorized distributors.

- Only use the INTRABEAM system within the framework of the intended use.
- Observe all the symbols and labeling on the packaging and on the components of the INTRABEAM System, see Page 43.

Required additional information

These Instructions for Use contain detailed explanations on the use of the INTRABEAM core system and information about the accessories. You can find additional explanations on the use of the INTRABEAM System in the Instructions for Use or product information, which are enclosed with the individual components of the INTRABEAM core system and the accessories.

- Also read the following Instructions for Use through completely before using the INTRABEAM core system:
 - INTRABEAM cart
Order Number, ZEISS, 000000-1501-288,
ZEISS Instructions for Use G-30-1704
 - NC32 INTRABEAM Floor Stand,
Order Number, ZEISS, 1096-945
ZEISS instructions for Use G-30-1476
 - INTRABEAM drape
Order Number, ZEISS, 32090-0000-000,
ZEISS instructions for Use G-30-1526
 - Radiation Shield, flat
Order Number, ZEISS, 304534-6000-595,
ZEISS instructions for Use G-30-1880
 - INTRABEAM spherical applicators
Order Number, ZEISS, 304534-6000-348,
ZEISS instructions for Use G-30-1478
 - INTRABEAM Needle Applicator,
Order Number, ZEISS, 304534-3001-200,
ZEISS instructions for Use G-30-1723

- INTRABEAM Surface Applicator Set 1.0 cm, 2.0 cm, 3.0 cm and 4.0 cm
Order Numbers, ZEISS, 304534-3001-530, -531, -532, -533;
ZEISS instructions for Use G-30-1858
- INTRABEAM Flat Applicator Set 1.0 cm, 2.0 cm, 3.0 cm, 4.0 cm,
5.0 cm and 6.0 cm;
Order numbers, ZEISS, 304534-3001-565, -566, -567, -568,
-569 and -570;
ZEISS instructions for Use G-30-1880
- UNIDOS E dosimeter from PTW, Type 10010,
Instructions for Use_en, PTW, 54513100_02
Instrucciones de manejo_es, PTW, 54519600_00
Gebrauchsanweisung_de, PTW, 54515100_02
- INTRABEAM Water Phantom,
Order Number, ZEISS, 304534-7000-600,
ZEISS instructions for Use G-30-1767
- PRS 500 user terminal
Instructions for Use from PENTA
Medical ETX-CD Zeiss manual English R 1.1 (on CD-ROM)

Electrical safety

- Switch off the device every time you are going to unplug it from the electricity network or connect it, have not used for an extended period of time, or want to clean its surface.
- Switch the device off at the power switch or unplug the power lead if you notice smoke, sparks or unusual smells. In this case, do not use the device again before it has been repaired by our Service department.
- If necessary, include the device in the safety measure "potential equalization". The potential equalization bolt (see chapter "System overview") can be used for connecting other active devices to the same ground potential or for redundant grounding to protective earth.
- Do not use any multiple sockets!
- Never open the device! The device contains freely accessible live components. If you remove the housing, you run the risk of electric shock.
- Do not touch the system if your body is electrostatically charged and the system is not grounded.
- Connect the device to the "emergency backup power supply" according to the country-specific regulations or guidelines.

- Transport & service*
- Over longer distances (e.g. removal, return for repair, etc.) the system must always be transported in the original packaging or in special return packaging. For details, please contact your dealer or the Carl Zeiss service team.

ZEISS Service This system is a high-grade technological product. We recommend that it is checked at regular intervals to ensure that it operates perfectly and safely.

- Have your device checked at regular intervals by Zeiss Service.

Carl Zeiss Meditec AG has appropriate Service Agreements available.

- In order to prevent any impairment of the device's safety due to age, wear, etc., the user must ensure that the device is subjected to the necessary safety checks (see section "Care and Maintenance").
- Label the device as non-functioning and notify our Service department if an error occurs that you cannot rectify using the table "Help with Troubleshooting".

Information on Zeiss Service The Service Manual "SM-30-8023 INTRABEAM System" is available to ZEISS Service for its work. This Service Manual has only been produced for Zeiss Service for their internal use.

Service Identification Program (SIP) If your unit or system has been included in the ZEISS "Service Identification Program" (SIP) it has been marked with a SIP label with a SIP number, see Page 46. The SIP number is a unique identification number assigned to your device or system. A product file is maintained by ZEISS for this SIP number.

- For any service requirements, please contact your local contact at the ZEISS subsidiary in your country. You will find their phone number, FAX number and email address on the SIP label.

Prerequisites for operation Our service representative or an expert authorized by us will install the system. Please ensure that the following requirements for operation are maintained:

- All mechanical connections (details in the user's manual) which are relevant to safety are properly connected and screw connections tightened.
- All cables and plugs are in good working condition.
- The voltage setting on the instrument conforms to the rated voltage of the line supply on site.
- The instrument is plugged into a power outlet which has a properly connected protective ground contact.
- The power cord being used is the one designed for use with this instrument.

- Do not force cable connections. If the male and female parts do not readily connect, make sure that they are appropriate for one another. If any of the connectors are damaged, have our service representative repair them.
- Do not pull on the cables. The plug connections are self-locking. Slide back the sleeve to separate the connection.

Before every use and after every system change

- Ensure that all specified "requirements for operation" are met.
- Go through the checklists, see Page 178.
- Note all symbols, notes and labels on the system.

Safety during the treatment

- Reduce the beam dose and the radiation period to what is absolutely necessary. Any type of radiation has a damaging effect on biological tissue.
- If an XRS x-ray source is mounted on the adapter plate of the INTRABEAM stand and an applicator is inserted into the adapter shaft of the adapter plate, this activates the visual safety barrier (interlock) and the INTRABEAM system can generate x-ray beams.
- Suitable radiation shielding must be used to protect patient tissue outside the target area and OR personnel during use.
- Only work with the system after it has been properly calibrated. For notes on calibration see Page 35.
- Never operate the system unattended.
- Never cover any ventilation openings.

After every use of the system

- Use the main power switch to switch off the system. The main power switch must always be off when the system is not in use.

Use of the INTRABEAM Core System

The operator is required to follow personnel selection, safety, training, verification, sterilization, and other requirements necessary to assure safety of both patients and operating personnel.

The XRS X-ray sources are calibrated in the ZEISS manufacturing plant. ZEISS supplies these sources with a calibration file, see Page 120. Without the associated calibration file, the XRS X-ray source cannot be operated in the INTRABEAM Core System.

- As a prerequisite for the use of the INTRABEAM Core System, the output of the XRS X-ray source will need to be calibrated by the manufacturer, Carl Zeiss Meditec AG, using a known reference which is traceable to a national standard.
- Switch the XRS X-ray emission source on only to treat a patient or when introducing it into one of the components used to verify the XRS X-ray emission source (PAICH and PDA).
- Only place the XRS X-ray radiation source into its mounting or attach an applicator if the X-ray radiation source is intended to emit X-radiation. Make sure you use suitable radiation shielding to protect untargeted patient tissue and operating personnel during use.
- Handle the IC ionization chamber carefully and do not touch the sensor area. The IC ionization chamber supplied with the INTRABEAM Core System is a sensitive and delicate device.
- To safeguard the electrical safety of the INTRABEAM Core System, always connect the system with the protective grounding in the building when it is in use.
- Do not remove any components of the INTRABEAM Core System as this may lead to malfunctions.
- Always wear protective clothing such as gloves, an apron or gown and an eye shield when handling contaminated components belonging to the INTRABEAM Core System.
- The INTRABEAM Core System is designed to be used by qualified personnel properly trained in its use. Please observe the instructions in the section "Requirements for the user", Page 85.
- You may only connect devices approved by ZEISS to the interfaces.
- You must turn the probe centering unit / PAICH Ionization Chamber holder PAICH completely ($\geq 360^\circ$) around the part that is to be measured, in order to be able to measure the out-of-center error accurately.

- The application software calculates the overall exposure time using the dose rate and the penetration depth. Manual entry of the treatment time is not possible.
- Please observe the instructions in the section "Handling, Transportation and Shipping, Page 239.
- Please observe the instructions in the section "Requirements for System Operation", Page 85.

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The owner/operator or his IT representative and every user must ensure the observance of the corresponding laws and data protection regulations for the respective state. The required technical and organizational prerequisites for this must be created so that

- the data in the system is protected against misuse, i.e. that any personal, protected data cannot be collected, processed or used without authorization.
- personal data is protected, i.e. that any individual details about personal or factual conditions of a certain or certain natural entities, i.e. all data which can be directly assigned to a person (employee, customer, supplier), such as marital status, type of occupation, religion, income, etc. cannot be passed on if it directly relates to a person.
- the data is protected during data processing, i.e. data is protected when storing (collection, recording or archiving), transferring (forwarding to parties outside the company), modifying (changing the content, including anonymization and aliasing), deleting and blocking (labeling to restrict further processing or use).
- data is also protected when used in another manner (e.g. internal dissemination).
- data forwarded to third parties outside the responsible center is protected by the recipient there.

Archiving of data

- The INTRABEAM System can only save data; it is not intended for the archiving of data.
- Data can be deleted by every subsequent user!
- In the event of repairs, all data is deleted in the user terminal. This means, all the data that was not archived previously, will be lost. For example, you can copy the data onto CD.

All users are responsible for archiving their own data.

If sufficient storage space is no longer available, the system informs you that files should be archived or deleted, if no longer required.

Safety devices of the INTRABEAM Core System

Optical interlock system

An optical interlock system has been incorporated into the XRS X-ray source to detect the type of radiation shielding integrated in the device into which the XRS is mounted.

If the XRS x-ray source is put into its holder and a Carl Zeiss applicator is connected, the optical interlock system closes and allows the generating of x-ray beams.

If x-ray beams are generated, an acoustic signal will be heard in the loud-speaker (6) of the control unit and the external warning lamp will light up.

The status of the interlock system is permanently displayed in the screen masks of the application software, see "Interlock Status" display field. The four possible interlock states are shown in the table below.

With the exception of the holder for the XRS x-ray source and the applicators, the generating of x-ray beams is only permitted, if devices with suitable radiation shielding are mounted on the x-ray source.

Interlock status	Generating of X-rays	Device
None	deactivated	None
non-shielded	deactivated	PAICH
shielded	permissible	PAICH with installed Ionization chamber (IC) PDA
treatment	permissible	XRS x-ray source and the applicator are properly mounted on the INTRABEAM stand

EMERGENCY OFF, control lamps and external interlock

- The control unit has an emergency OFF switch.
- The control unit has a loudspeaker.
 - When all interlock connections have been closed (the system is in the "Ready" state) long beeps are heard.
 - During X-radiation emission, short beeps are heard.
- The XRS X-ray source is equipped with an indicator light which lights up orange when X-radiation is being generated.
- The control unit visually indicates the production of X-radiation. The control console has indicator lamps for the following:
 - Power On, Color: white
 - X-rays On, Color: Orange
 - X-ray Ready, Color: Green
 - X-ray Pause, Color: Blue

The indicator lamps are described on Page 54.

- There is an external interlock connection on the reverse side of the control unit which can be used both as an external interlock connection as well as a remotely controlled display of X-ray generation.

For explanations of the external interlock connection and the X-ray display, see Page 250.

The INTRABEAM Core System is designed to protect both the patient and the user in the event of electrical problems by diverting currents over external surfaces which are connected to the earth of the building electrical system via the power supply cord.

A potential equalization connection can be found on the rear side of the control console. This connection can be used for connecting other active devices at the same ground potential as the INTRABEAM Core System.

The INTRABEAM Core System is equipped with electrical fuses to reduce possible hazards due to an electrical system failure. In the event that a fuse needs to be replaced, only replace it with a fuse that has the ratings indicated on the rear panel of the control console.

Calibrations

The use of the INTRABEAM Core System requires that the user arranges for the following components to be calibrated:

- XRS X-ray radiation source
- Ionization chamber (IC)

ZEISS recommends annual calibration of these components.

PTW recommends calibration of the UNIDOS E dosimeter every two years.

- Follow the additional instructions in the section "Calibration rules", Page 209.

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When the XRS X-ray source is used with the PAICH (Probe Adjuster/Ionization Chamber Holder) to straighten the probe, the high voltage for the XRS is not turned on and no x-rays are generated.

When the XRS X-ray source is used with the PAICH (Probe Adjuster/Ionization Chamber Holder) to check the X-ray output with the IC (ionization chamber) installed in the PAICH, the XRS generates X-rays.

Due to the IC shielding, the PAICH provides adequate radiation shielding when used as described in this manual.

When the XRS X-ray source is used with the PDA photodiode array to check isotropy, the XRS generates X-rays.

The PDA photodiode array provides adequate radiation shielding when used as described in this manual.

Safety during treatment

Physician using the system should always carefully weigh up, whether radiation with the INTRABEAM core system is suitable when selecting patients for radiation with the INTRABEAM core system, as patients who have severe edema, bleeding or severe calcification in the treatment area.

X-rays are only produced if the XRS x-ray source is activated.

The control unit monitors the radiation output during the treatment with the internal IRM x-ray monitor. Upon reaching the prescribed radiation dose, the internal IRM x-ray monitor stops the treatment.

The IRM internal radiation monitor detects the radiation that the patient is exposed to during treatment.

At the beginning of the treatment, the application software calculates the total exposure time based on the entered dose and the penetration depth. If the internal x-ray monitor does not stop the treatment (for example, because its counting rate has dropped), the treatment time will be used as a discontinuation criterion: If the initially calculated total exposure time is exceeded by 10%, the timer function stops the treatment.

- It is necessary to use a sterile applicator, when you introduce the XRS x-ray source probe into a patient.
- To ensure that the treating staff is not exposed to any excessive radiation dose, the ambient radiation should be monitored during the treatment with a suitable test device.
- Wearable radiation shielding should be used to protect the operating personnel.
- Switch off the x-ray beam before you remove the probe on the XRS x-ray source out of the patient.

An unavoidable side effect during the production of x-ray radiation with the probe of the XRS s-ray source is the simultaneous generating of thermal energy and its dissipation into the surrounding tissue of the patient. The surface temperature of the spherical applicators could exceed 41°C during the treatment session.

Electrical features of the INTRABEAM Core System

The control console is equipped with various electrical connections on the back of the unit. Only devices authorized by ZEISS may be connected here.

The system is equipped with electrical fuses to reduce possible hazards due to an electrical system failure and to ensure sustained protection against fire.

In the event that a fuse needs to be replaced, only replace it with a fuse of the same type and with the same ratings. Observe the data on the rear panel of the relevant device.

Never use a power cord, a power supply with isolating transformer or a power cord adapter that is not compliant with IEC 60601-1 and does not maintain earth grounding of the INTRABEAM Core System.

To ensure electrical safety, the INTRABEAM Core System must always be grounded to the protective earth of the room when it is in use. Always use the power cord that is supplied with the INTRABEAM Core System, and ensure proper grounding of the mains supply.


ZEISS recommends an additional ground connection for the INTRABEAM Core System via the potential equalization of the control unit, the INTRABEAM floor stand or the isolating transformer in the equipment trolley.

The maximum voltage supplied to the XRS X-ray source is ± 15 V.

Any hazard caused by exposure to high voltage is prevented by generating the high voltage in a well-insulated housing, i.e. the electronic system is surrounded by the metal case of the XRS X-ray source. The metal cladding on the XRS X-ray radiation source is designed such that it can be connected to the earth in the room via the XRS X-ray radiation source cable and the control console power cable.

- Take appropriate precautions when using other grounded or ungrounded devices in conjunction with the INTRABEAM Core System.

A potential equalization connector (1) is located on the rear panel of the PRS 500 control console. The connection of the PRS 500 control console to the facility power supply ensures that the potential equalization connector is connected via the power cord to the protective earth of the control console and the protective earth conductor of the facility electrical system.

You can use the potential equalization connector for placing other active devices at the same ground potential as the INTRABEAM Core System. Electrical potentials between the INTRABEAM Core System and other grounded devices in the vicinity can be minimized by connecting the potential equalization connectors of all other devices to the potential equalization connector on the PRS control console.

The power inlet on the PRS control console is provided with a stirrup (2) to prevent inadvertent unplugging of the power plug.

- Always fold up stirrup (2) to support the power plug when using the PRS 500 control console.

The INTRABEAM Core System components are not defibrillation-proof. Therefore remove the XRS X-ray source from the patient before defibrillating.

Fire, anesthetics and oxygen

Inflammable substances such as anesthetics and oxygen should not be used near activated electrical components of the INTRABEAM core system (control unit, user terminal).

Do not switch on components of the INTRABEAM Core System powered by electricity while they are in the vicinity (within 25 cm) of flammable anesthetics or an atmosphere enriched with oxygen or nitrogen monoxide.

EMC (Electromagnetic compatibility)

The INTRABEAM System meets the EMC requirements of IEC 60601-1-2. You must not use the INTRABEAM System alongside other ME devices and systems. If use in the immediate vicinity of other devices cannot be avoided, you must constantly monitor the INTRABEAM system during operation to ensure compliance with the intended use.

- When operating the INTRABEAM system, please observe the following precautions regarding EMC:
 - Only use spare parts approved by ZEISS for this device.
 - Do not use any portable or mobile RF communication equipment in the vicinity of the device as this may impair the device's function.
 - Do not use any mobile phones near the device. They represent a potential risk of causing the incorrect operation of medical devices.
 - Malfunctions may occur, depending on a variety of local factors. They cannot be predicted and can by no means be estimated.
- Observe the EMC guidelines in the Section "EMC - Electromagnetic Compatibility", see Page 252ff.

Cleaning, disinfection, sterilization

- To protect patients and operating personnel against infections, only use components during the treatment that have been carefully cleaned, disinfected and sterilized.
- Use pre-sterilized, disposable drapes.

Sterility cannot be guaranteed or damage can occur if non-permissible methods for cleaning, disinfecting and sterilizing components are used or if components are cleaned, disinfected and sterilized that are not intended for this. This may lead to excess radiation.

- Only use the permissible cleaning, disinfection and sterilization methods described here.
- Do not clean, disinfect and sterilize any disposable products.

Further details and information on permissible cleaning, disinfecting and disinfection methods can be found on Page 234 ff.

INTRABEAM Components and Accessories	Washing	Rinsing	Disinfection	Sterile single-use product or reconditionable product
NC32 INTRABEAM Floor Stand	X	-	X	-
INTRABEAM cart	X	-	X	-
PRS 500 control unit	X	-	X	-
User terminal	X	-	X	-
Keyboard	X	X	X	-
XRS 4 X-ray source	X	-	X	-
PDA photo diode array	X	-	X	-
PAICH	X	-	X	-
UNIDOS E dosimeter*	X	-	X	-
Ionization chamber IC*	X	-	X	-
XRS, PDA/PAICH cables	X	X	X	-
V-block	X	X	X	-
X-block	X	X	X	-
INTRABEAM Water Phantom	X	-	X	-
Transport and storage container	X	X	X	-
INTRABEAM shields	-	-	-	Sterile Single-use product
INTRABEAM spherical applicators	X	X	X	steam-sterilizable
INTRABEAM Surface Applicator set	X	X	X	steam-sterilizable
INTRABEAM Flat applicator set	X	X	X	steam-sterilizable
INTRABEAM needle applicator	-	-	-	Sterile Single-use product
INTRABEAM drapes	-	-	-	Sterile Single-use product

* If you have to clean and disinfect the UNIDOS E Dosimeter and the IC Ioni-

zation Chamber, please refer to the corresponding Instructions for Use by PTW, Freiburg.

Disposal of consumables (single-use material)

Used consumables, such as drapes, shields and needle applicators are infectious waste, which must be disposed of properly.

- Follow the relevant directives of your hospital.
- Observe the relevant national codes and regulations.

Disposal of components of the INTRABEAM Core System

Used or obsolete INTRABEAM Core System components can be returned to the manufacturer ZEISS after cleaning.

- For the disposal of such components, please get in touch with your local contact at the ZEISS subsidiary in your country.


1. Introduction
 2. Background
 3. Methodology
 4. Results
 5. Discussion
 6. Conclusion
 7. References
 8. Appendix
 9. Index
 10. Table of Contents
 11. Figure 1
 12. Figure 2
 13. Figure 3
 14. Figure 4
 15. Figure 5
 16. Figure 6
 17. Figure 7
 18. Figure 8
 19. Figure 9
 20. Figure 10
 21. Figure 11
 22. Figure 12
 23. Figure 13
 24. Figure 14
 25. Figure 15
 26. Figure 16
 27. Figure 17
 28. Figure 18
 29. Figure 19
 30. Figure 20
 31. Figure 21
 32. Figure 22
 33. Figure 23
 34. Figure 24
 35. Figure 25
 36. Figure 26
 37. Figure 27
 38. Figure 28
 39. Figure 29
 40. Figure 30
 41. Figure 31
 42. Figure 32
 43. Figure 33
 44. Figure 34
 45. Figure 35
 46. Figure 36
 47. Figure 37
 48. Figure 38
 49. Figure 39
 50. Figure 40
 51. Figure 41
 52. Figure 42
 53. Figure 43
 54. Figure 44
 55. Figure 45
 56. Figure 46
 57. Figure 47
 58. Figure 48
 59. Figure 49
 60. Figure 50
 61. Figure 51
 62. Figure 52
 63. Figure 53
 64. Figure 54
 65. Figure 55
 66. Figure 56
 67. Figure 57
 68. Figure 58
 69. Figure 59
 70. Figure 60
 71. Figure 61
 72. Figure 62
 73. Figure 63
 74. Figure 64
 75. Figure 65
 76. Figure 66
 77. Figure 67
 78. Figure 68
 79. Figure 69
 80. Figure 70
 81. Figure 71
 82. Figure 72
 83. Figure 73
 84. Figure 74
 85. Figure 75
 86. Figure 76
 87. Figure 77
 88. Figure 78
 89. Figure 79
 90. Figure 80
 91. Figure 81
 92. Figure 82
 93. Figure 83
 94. Figure 84
 95. Figure 85
 96. Figure 86
 97. Figure 87
 98. Figure 88
 99. Figure 89
 100. Figure 90
 101. Figure 91
 102. Figure 92
 103. Figure 93
 104. Figure 94
 105. Figure 95
 106. Figure 96
 107. Figure 97
 108. Figure 98
 109. Figure 99
 110. Figure 100
 111. Figure 101
 112. Figure 102
 113. Figure 103
 114. Figure 104
 115. Figure 105
 116. Figure 106
 117. Figure 107
 118. Figure 108
 119. Figure 109
 120. Figure 110
 121. Figure 111
 122. Figure 112
 123. Figure 113
 124. Figure 114
 125. Figure 115
 126. Figure 116
 127. Figure 117
 128. Figure 118
 129. Figure 119
 130. Figure 120
 131. Figure 121
 132. Figure 122
 133. Figure 123
 134. Figure 124
 135. Figure 125
 136. Figure 126
 137. Figure 127
 138. Figure 128
 139. Figure 129
 140. Figure 130
 141. Figure 131
 142. Figure 132
 143. Figure 133
 144. Figure 134
 145. Figure 135
 146. Figure 136
 147. Figure 137
 148. Figure 138
 149. Figure 139
 150. Figure 140
 151. Figure 141
 152. Figure 142
 153. Figure 143
 154. Figure 144
 155. Figure 145
 156. Figure 146
 157. Figure 147
 158. Figure 148
 159. Figure 149
 160. Figure 150
 161. Figure 151
 162. Figure 152
 163. Figure 153
 164. Figure 154
 165. Figure 155
 166. Figure 156
 167. Figure 157
 168. Figure 158
 169. Figure 159
 170. Figure 160
 171. Figure 161
 172. Figure 162
 173. Figure 163
 174. Figure 164
 175. Figure 165
 176. Figure 166
 177. Figure 167
 178. Figure 168
 179. Figure 169
 180. Figure 170
 181. Figure 171
 182. Figure 172
 183. Figure 173
 184. Figure 174
 185. Figure 175
 186. Figure 176
 187. Figure 177
 188. Figure 178
 189. Figure 179
 190. Figure 180
 191. Figure 181
 192. Figure 182
 193. Figure 183
 194. Figure 184
 195. Figure 185
 196. Figure 186
 197. Figure 187
 198. Figure 188
 199. Figure 189
 200. Figure 190
 201. Figure 191
 202. Figure 192
 203. Figure 193
 204. Figure 194
 205. Figure 195
 206. Figure 196
 207. Figure 197
 208. Figure 198
 209. Figure 199
 210. Figure 200
 211. Figure 201
 212. Figure 202
 213. Figure 203
 214. Figure 204
 215. Figure 205
 216. Figure 206
 217. Figure 207
 218

Warning and information labels

See the warning and information labels!

If you notice that any label is missing on your device or has become illegible, please contact us or one of our authorized representatives. We will provide you with replacement.

PRS 500 Control Console rating plate


The rating label provides the following information:

- Symbol for "Manufacturer" 
- Manufacturer (company name) Carl Zeiss Meditec AG
- Manufacturer's address Goeschwitzer Strasse 51-52
07745 Jena, Germany
- Serial number symbol  640210xxxx
- IEC number IEC 60601-2-8:1997
- Device name Control Unit PRS500
- IEC number IEC 60601-2-8:2010
- Reference number symbol  304534-7000-338
- Product classification P30600001
- Rated voltage 100 V - 240 V
- Line frequency range 50 Hz - 60 Hz
- Connected load 60 VA
- Symbol for ionizing radiation 
- CSA label 
- CE labeling with number of notifying body  0297
- Further technical data about the PRS 500 control unit can be found in the section "Technical Data" on Page 248.
- The technical data for the user terminal is located in the section "Technical Data" on Page 249.

Location of the rating label

The rating label (1) is attached to the rear panel of the PRS 500 control console.

Symbols and labeling on the back of the control unit PRS 500

Potential equalization connector
another electrical device can be connected to the same mass potential as the INTRABEAM core system here.

ON (power supply, connection with line power)

OFF (power supply, no connection with line power)

REPLACE ONLY WITH TYPE T1.25A E/250V FUSES
(Only replace fuse with type T1.25A E/250V)

Explanation about fuse type "T1.25A E/250V":

- time delay type, slow-blow
- rated current, 1.25 A
- breaking capacity, increased
- Rated voltage, 250 V

1.25A

250V

Notice "Observe Instructions for Use"


Observe Instructions for Use or accompanying documents

Device classification, User part Type B pursuant to EC 60601-1


Warning label

The WARNING label contains the following information:

- 1 Grounding reliability can only be achieved when the equipment is connected to the receptacle marked "hospital only" or "hospital grade".
(Reliable earthing is only ensured if the device is connected to the outlet labeled "hospital only" or "hospital grade".)
- 2 This graphic symbol means:
"Caution, see accompanying documentation"
- 3 Symbol for Ionizing radiation
- 4 This X-ray unit may be dangerous to patient and operator unless safe exposure factors and operating instruction are observed.
(In order that the x-ray unit does not pose a danger for patients and operating personnel, it is mandatory the conditions for a safe delivery of radiation and the Instructions for Use be observed.)
- 5 Location for "WARNING" signs (see Figure)

Symbols and inscriptions laterally and at the bottom


The SIP label provides the following information:

- Manufacturer (company name)
- Manufacturer's contact information
i.e. phone number, fax number and email address of the local contact of the national ZEISS sales organization.
- SIP number
A unique identification number assigned to your device or system.
A product file is maintained by Carl Zeiss for this SIP number.
- Location of the SIP label (1)
The SIP label is attached to the right side of the PRS 500 control console.


Sign "Observe disposal regulations" (2)

This symbol means that electrical or electronic devices must not be disposed of as normal domestic waste.


"Date of manufacture" plate (3)

This label indicates the date of manufacture (year - month) of the device.

Label on the XRS 4 X-ray source

Key to symbols

This symbol indicates the emission of X-ray radiation.

On radiological devices
this symbol indicates the X-ray source module.

"XRS measured" label

This label is used to mark the XRS X-ray source which was measured with an INTRABEAM Flat Applicator and an INTRABEAM Surface Applicator. A measurement report is issued to the customer.


Rating label of the PRS 500 user terminal

The rating label provides the following information:

- | | |
|-------------------------------|--|
| - Manufacturer (company name) | Carl Zeiss Meditec AG |
| - Symbol for "Manufacturer" | |
| - Manufacturer's address | Goeschitzer Strasse 51-52
07745 Jena, Germany |
| - Serial number symbol | |
| - Product classification | P30600001 |
| - Device name | PRS 500 Terminal (Intrabeam) |
| - Rated voltage | 100 V - 240 V |
| - Connected load | max. 100 VA |
| - Line frequency range | 50Hz - 60Hz |
| - CSA label | |
| - CE label | |

Location (1) of the rating label

System overview

System description	50
Main components of the INTRABEAM Core System	50
Verification components.....	50
Components for the protection of the INTRABEAM Core System	51
Accessories of the INTRABEAM system	51
 Principal components	 53
PRS 500 control console	53
User terminal.....	58
XRS 4 X-ray radiation source and IRM internal X-ray radiation monitor	64
 Protection components	 66
X-block.....	66
V-block guide	67
Probe protector	68
 Verification devices	 70
Probe adjuster/ionization chamber holder PAICH	70
PDA photodiode array	72
Ionization chamber (IC).....	72
Cables - color coding	74
UNIDOS E dosimeter.....	75
Water phantom, optional verification accessories	76
Transport and storage container.....	78
 Accessories for the INTRABEAM core system	 80
INTRABEAM cart	80
NC32 INTRABEAM Floor Stand	81

System description

The INTRABEAM PRS 500 Core System with XRS 4 comprises several components.

Different system configurations are required for different types of treatment.

Main components of the INTRABEAM Core System

- PRS 500 control console
- XRS 4 X-ray source with internal radiation monitor (IRM)
- User terminal for configuration and operation of the INTRABEAM Core System

For illustrations see Page 13.

Verification components

- V-block guide
- Photodiode array (PDA)
- Probe adjuster/ionization chamber holder (PAICH)
- Ionization chamber with shielding
- UNIDOS E dosimeter

For illustrations see Page 14.

Accessories for checking the system

INTRABEAM Water Phantom

Information for INTRABEAM Water Phantom, see the corresponding User guide, see list on Page 26.

Components for the protection of the INTRABEAM Core System

- Probe protector on XRS
For illustration see Page 68.
- X-block, for illustration see Page 66.

Accessories of the INTRABEAM system

Since the products are not approved in all countries, please contact with respect to the approval of the local partners.

Information on the accessory components of the INTRABEAM system, refer to the corresponding User instructions. Overview, see Page 26.

- INTRABEAM Spherical Applicators
- INTRABEAM Needle Applicator
- INTRABEAM Flat Applicator Set
- INTRABEAM Surface Applicator Set
- INTRABEAM Radiation Shield, Flat
- INTRABEAM Drape
- NC32 INTRABEAM Floor Stand
- INTRABEAM cart

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Principal components

PRS 500 control console

The control console receives job packets from the user terminal for verification and during treatment.

The control console executes these job packets and independently monitors the process. It also displays the current data during this process.

The control console supplies raw data to the user terminal which uses this data to generate treatment and verification records.

The control console consists of the following major components:

- Software and electronics to monitor and control various system components
- Display for setting parameters, indicator lights and buttons
- Power supply unit that supplies the low voltage required by the XRS
- Air pressure sensor

Controls, displays, connections

Front of the PRS 500 control console

WARNING:

If you do not precisely know the functions of the control and display elements, you will not be able to properly operate the INTRABEAM system. Faulty operation may lead to radiation damage.

- Make yourself familiar with the functions of the controls and displays.

1 Emergency stop button

This large red button can be used at any time to switch off the control console if a critical situation develops.

If you press the emergency off button in a critical situation during treatment in which X-ray emission is in progress and the orange indicator light "X-RAYS ON" is illuminated, treatment is stopped immediately and X-ray emission also stops. The treatment parameters are saved before treatment stops.

If you abort treatment by pressing the Emergency Stop button, unless you continue to treat the same patient by using the "Continue with previous patient" option, you must re-verify the XRS X-ray source. In that case, the remaining dose can be administered within 90 minutes without the need to re-verify the XRS x-ray source, see the "Enter prescription" screen Page 167.

To switch the PRS 500 control console on again, you must first unlock the Emergency Stop button (turn the button and pull it outwards).

2 POWER ON indicator light (color: white)

Illuminated when the control console power is on.

3 Control console display

Displays data required to monitor the INTRABEAM Core System components. The contents of the control panel display vary according to the current procedure, see items 8 - 15.

4 X-RAYS ON indicator light (color: orange).

As on the XRS x-ray source.

Illuminated when the XRS x-ray source is generating X-rays

5 READY indicator light (color: green)

Illuminated when all preparatory activities have been completed and the system is ready to generate X-rays.

6 PAUSE indicator light (color: blue)

Illuminated when X-ray generation is temporarily disabled and the system is in the PAUSE mode.


7 "PAUSE" button

By pressing the "PAUSE" button you temporarily interrupt X-ray emission. The control parameters are retained. The PAUSE indicator light is illuminated.



Press the "Resume" button in the treatment screen on the user terminal (application software) to resume X-ray emission.

8 - 15 Contents of the control panel display vary according to current procedure:

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Control panel display during system verification (SQA):


Rear panel of the PRS 500 control console

- 1 Connection of an external interlock and external warning light.
To connect the cable for an external interlock switch (eg door contact), please use the enclosed "Interlock Control Plug".
The interlock requirements must be satisfied in order to produce X-rays. During the emission of X-rays, an audio signal is heard from loudspeaker (6) and the external warning light is illuminated.

The system is delivered with a dummy plug (2) inserted in the connector for the external interlock cable as this is an optional connection possibility. If the devices connected to the external interlock connector are also connected to line voltage, these devices must meet the electrical safety requirements specified in IEC 60601-1.
- 2 Dummy plug for the external interlock connector
- 3 Potential equalization connector
Use the potential equalization bolt to incorporate the PRS 500 control console into the potential equalization measures. You can use the potential equalization bolt to connect separate pieces of equipment electrically, placing them at the same ground potential.
- 4 Fuse
The PRS 500 control console's internal circuitry is protected against power surges via a fuse. Replacement of the fuse is described on Page 233.

Observe the note on the instrument label plate:
"REPLACE ONLY WITH TYPE T1.25A E/250V FUSES"
- 5 Power switch
In the "I" position, the line voltage is switched on for the PRS 500 control console. In the "O" position, the line voltage is switched off for the PRS 500 control console.
- 6 Loudspeaker
An audio signal is heard from the loudspeaker during the emission of X-rays.

- 7 Power inlet socket, connector for IEC power plugs
Permissible line voltages 100 - 240 VAC, 50 - 60 Hz
The stirrup on the power inlet socket prevents inadvertent unplugging of the power plug.
- 8 Rating label
See also the "Warning labels and notes" chapter on Page 43.
- 9 Serial interface port "COM 1"
This is the RS232 interface port for the user terminal.
- 10 PDA/PAICH connector
PDA/PAICH cable connector, yellow color code.
- 11 Only for predecessor systems:
ERM connector
ERM cable connector, blue color code.
- 12 XRS connector
XRS cable connector, green color code.

User terminal

PRS 500 user terminal

The user terminal is used to configure and operate the INTRABEAM Core System.

The user terminal comprises:

- Computer, integrated in the monitor housing
- Touchscreen monitor to activate functions
- Membrane keyboard with touchpad (1), to enter figures and text
- Drives for floppy disks and DVD/CD-ROM
- Various ports, including two RS232 interfaces.

For full-mode operation (online mode) of the INTRABEAM Core System, you have to connect the user terminal to the control console and to the UNIDOS dosimeter, each one via one of the RS232 interfaces.

The application software installed on the computer enables you to:

- capture, activate, deactivate and reformat data of the components used
- check the system
- enter treatment parameters
- verify treatment parameter settings
- perform treatments
- review treatments
- document and save verification and treatment

1 "POWER" button

This button enables you to switch between the operating mode and standby mode of the user terminal.




2 Touchpad (cursor control), separate mouse buttons(available for selection)

3 USB port




You can only connect a USB stick. Other USB devices must not be connected for safety reasons.

4 Drives for floppy disks and DVD/CD-ROM

5 PCMCIA slot with external memory card for database backup

- 6 Power source unit.
consisting of main on/off switch, fuse holder and cooler jack
- 7 Connection bolt for potential equalization
- 8 RS232 interface - COM 2
- 9 RS232 interface - COM 1
- 10 Connection for keyboard with PS2 interface
- 11 Connection for keyboard with USB interface

Operating the PRS 500 user terminal


- 1 <Increase screen brightness> button
- 2 <Reduce screen brightness> button
- 3 <Deactivate touch function> button
The touch function can be deactivated to avoid unintended faulty operation, e.g. while cleaning the touchscreen. When the touch function has been deactivated, the "TOUCH" indicator lamp (4) is lit.
Press this button again to reactivate the touch function.
- 4 The "TOUCH" indicator lamp
illuminates when the touch function is deactivated.
- 5 The "HDD" indicator lamp
is the heartbeat indicator of the user terminal: it blinks while the operating system is running.
- 6 The "POWER" indicator lamp
illuminates when the user terminal is in the operating mode.
- 7 <POWER> button
This button enables you to switch between the operating mode and standby mode of the user terminal.

Normally, the user terminal is not separated from line power (no galvanic isolation).

To turn on the user terminal (switch it to the operating mode):
Press the <POWER> button for about 3 seconds until the "POWER" indicator lamp (6) illuminates. The operating system is started.

To turn off the user terminal (switch it to the standby mode):
Press the <POWER> button for about 3 seconds until the user terminal emits an acoustic signal.
Release the <POWER> button and briefly press it again. The "POWER" indicator lamp (6) goes out and the user terminal switches to the standby mode.

1

2

3

4

5

6

7

Operating the touchscreen

- Operate the touchscreen with your fingers.
- Never use any hard, sharp, pointed, wooden, metal or plastic objects such as pencils, ball-point pens or knives as this would destroy the touchscreen.

XRS 4 X-ray radiation source and IRM internal X-ray radiation monitor

High voltage is generated by the XRS X-ray radiation source from the low DC voltage supplied by the control console and is used to direct an electron beam into the X-ray probe.

The high voltage generator in the XRS X-ray source XRS is completely surrounded by a grounded metal casing to preclude the possibility of patient injury from the high voltage. To use the XRS X-ray radiation source it must be connected to the PRS control console with the XRS interface cable. Since the high voltage is generated by the XRS X-ray radiation source itself, there is no high-voltage cable leading to the XRS X-ray radiation source.

The electron beam strikes a target at the end of the XRS probe, generating X-ray radiation. The output of the XRS X-ray source is a near-spherical radiation field emitted from the tip of the probe.

The probe tip is made of beryllium (Be), a material which is permeable to X-radiation. The exterior surface of the probe tip is provided with a protective coating.

The XRS X-ray radiation sources must be calibrated on a regular basis. Have the XRS X-ray radiation source calibrated by ZEISS every year. If shorter calibration intervals are stipulated by national codes and regulations, these must be obeyed.

CAUTION

The probe tip of the XRS X-ray radiation source contains beryllium (Be). The exterior surface of the probe tip is provided with a protective coating. If the coating is damaged, skin and clothing may be contaminated.

- Avoid any direct skin contact with beryllium and immediately return the XRS X-ray radiation source to the manufacturer.
- First aid measures: Wash with lots of water after skin contact with beryllium. Remove contaminated clothing.


CAUTION

The probe tip of the XRS X-ray radiation source contains material incompatible with the human body. If the user inserts the XRS probe directly into the tissue, incompatible material comes into contact with the wound tissue. This may lead to allergic reactions.

- Never insert the XRS probe directly into the tissue.
- Always mount an applicator on the XRS probe.

- 1 Cable connection for the XRS 4 X-ray source
- 2 Shaft (deflector unit)
- 3 Probe tip (location where X-radiation is produced)
- 4 X-ray probe
- 5 X/Y marking for alignment of the PDA photodiode array
- 6 IRM internal X-ray radiation monitor

Internal radiation monitor (IRM)

The IRM is the primary monitor of the delivery of radiation to the patient.

The X-rays produced at the tip of the X-ray probe feature a spherical radiation characteristic. As a result, part of the radiation re-enters the XRS X-ray source through the probe. The IRM internal radiation monitor detects these X-rays in the XRS.


Before inserting the probe into the patient for treatment, check if the IRM detects any radiation output. The radiation output measured by the IRM is then used during treatment as an indirect measure of radiation delivered at the probe tip.

Refer to Page 114 for more information on calibrating the X-ray output of the XRS.

Protection components

X-block

The X-block is a mount for the XRS X-ray source, which holds the XRS in a stable position.


V-block guide

To determine whether the XRS X-ray source is functioning correctly, insert the XRS probe in the following verification devices:

- PAICH probe adjuster/ionization chamber holder
- PDA photodiode array

The V-block guide facilitates the insertion of the XRS X-ray source into the verification devices and prevents damage to the XRS probe.


Probe protector

The probe protector (1) is used to protect the XRS probe (4) from mechanical damage.

If you do not use the XRS X-ray source (4) after removing it from its storage case, you must slide the probe protector onto the XRS probe and mount it on the XRS using locking screw (2).

- 1 Probe protector
- 2 Probe protector locking screw
- 3 Mounting shaft for applicators
- 4 XRS X-ray source
- 5 INTRABEAM stand, mount for XRS X-ray source

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. A vertical red margin line runs down the right side of the page. The paper appears to be from a notebook or a standard composition book. There are no markings, text, or drawings on the page.

Verification devices

Probe adjuster/ionization chamber holder PAICH

The PAICH probe adjuster/ionization chamber holder has the following major functions:


- It measures the straightness of the XRS probe and, with predecessor systems, the dilators
- It enables the user to adjust (straighten) the XRS probe
- It provides a radiation-shielded location for placement of an ionization chamber during measurement of X-rays from the XRS.

The PAICH is sometimes referred to as the probe adjuster in reference to the first two functions.

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
The axis of the diode light source and photo detector is located so that the edge of the XRS probe partially interrupts the beam. If the XRS probe is straight, rotation of the PAICH will cause only minimal change in the light intensity transmitted to the photo detector. However, if the probe is not straight and the PAICH is rotated, the light signal intensity will vary.

The PAICH can then be positioned in such a way that the plunger can be used to straighten the probe.

The use of the PAICH as an ionization chamber holder is described in the IC (ionization chamber) section below.

1 Mount for the IC ionization chamber

This mount contains a temperature sensor. The INTRABEAM Core System requires the exact temperature to compute the dose correction during the "PAICH Output Check", see display window in the "PAICH Output Check" screen, see Page 157.

PDA photodiode array

The PDA photo diode array is used to measure and, if needed, adjust the isotropy of the radiation emission from the XRS probe tip.

The PDA photodiode array incorporates five photodiodes, placed in a centered position on the four side faces and the top face of a cube, such that all are equidistant from the center of radiation of the XRS probe tip.

The signals from the PDA are displayed in the screen and are used to measure the distribution (isotropy) of X-rays emitted from the tip of the probe. The isotropy of the radiation field is automatically adjusted on the basis of these measured values.

1 Connector (cable from control console)

2 Opening for insertion of the XRS probe

Ionization chamber (IC)

The IC ionization chamber supplied with the Core System is a sensitive and delicate device. Do not sterilize the IC ionization chamber.

WARNING

Touching the sensor window (1) can damage the IC ionization chamber to an extent that it is no longer operational. This may lead to over-irradiation.

- Do not touch the sensor window; always hold the IC ionization chamber at its side or cable.

The IC must be properly installed in the IC shield to ensure proper operation and radiation safety.

The IC supplied by Carl Zeiss has been calibrated in units of gray per coulomb (Gy/C) by a certified body. This calibration is traceable to the German National Metrology Institute (PTB) and optionally to the US National Institute for Standards and Technology (NIST). It is the responsibility of the operator to calibrate the IC to a standard specified by regulations that apply to the customer site. If no regulations apply, a standard of the operator's choice should be selected.


**WARNING**

If the operator does not observe the calibration intervals specified for the IC ionization chamber, the measurements performed during the output check can be incorrect. This may lead to incorrect irradiation.

- Have your IC calibrated at regular intervals.

The IC is used with a dosimeter to measure the dose rate from the XRS X-ray source. Prior to treatment, the IC is used for a verification process to check the output of the XRS 4 and to adjust the IRM. This corrects the calibration data and ensures the emission of an exact dose.

- 1 IC (ionization chamber) with detector window
- 2 IC shield
- 3 Connecting cable (to dosimeter)

Cables - color coding

The INTRABEAM Core System cables are color coded at the end connecting to the control console to match their proper receptacles.

They have gray connectors at the ends that connect to the respective components.

Color coding		
Cable	Connection to the components	Connection to the PRS 500 control console
XRS	gray	green
PDA/PAICH	gray	yellow

UNIDOS E dosimeter

The INTRABEAM Core System comes equipped with a PTW dosimeter (UNIDOS E model). The dosimeter is used to measure the current generated by the IC when the IC is exposed to X-rays. The UNIDOS E dosimeter is suitable for continuous operation.


Used with an XRS that is emitting X-rays, the IC generates typical currents of 1 to 60 picoamperes, when biased with +300 V DC. (This bias is automatically adjusted by the INTRABEAM PRS 500 system).

For instructions on using the UNIDOS E dosimeter from PTW, refer to the relevant manual. Read the entire PTW UNIDOS E dosimeter manual before using the dosimeter.


After powering the system, allow the UNIDOS E dosimeter to warm up for at least 5 minutes before you start the measurements.

Water phantom, optional verification accessories

The INTRABEAM Water Phantom is used to measure the depth dose curve and isotropy of the INTRABEAM system's XRS X-ray source in order to check if the measured values lie within the expected range.

- Carefully read the respective user manual before using the water phantom.

(Blank page, for your notes . . .)

[illegible]

Transport and storage container

Special transport and storage cases are available for the following INTRA-BEAM Core System components:

- XRS X-ray radiation source , probe protector and cable
- Special centering device / PAICH, PDA and cable.

The transport and storage cases can be used for storage and local transport in the hospital.

NOTE!

If the components are not stored in the relevant transport or storage case, they can be damaged and contaminated and therefore become unusable.

- Always store the components in their transport and sterilization case when you are not using them.

NOTE!

If the components of the INTRABEAM system are incorrectly inserted in their transport or storage case, they can be damaged and contaminated during in-house transport and therefore become unusable.

- Make sure you know how to fill the transport and storage cases correctly.


NOTE!

The transport and storage containers are not suitable for mailing the components. This may damage the components and make them unusable.

- Ask the manufacturer ZEISS for suitable shipping packaging.

Special sterile containers are also provided for spherical applicators. Please see the Instructions for Use for spherical applicators.

Transport and storage case for XRS, probe protector and cable**Transport and storage case for PAICH, PDA and cable**

Accessories for the INTRABEAM core system

INTRABEAM cart


The INTRABEAM cart is used exclusively for transport and storage of the INTRABEAM system components. It serves as a system rack and control station for INTRABEAM system verification and for the treatment process.

- Before using the INTRABEAM cart, read the user manual for the isolating transformer and, if necessary, the user manual of the isolation monitor.

Further information on the INTRABEAM cart


- How to relocate the INTRABEAM cart and how to correctly position it is described on Page 87.
- Cleaning and disinfecting the INTRABEAM cart is described on Page 238.

NC32 INTRABEAM Floor Stand

Detailed information on this accessory of INTRABEAM Core System, refer to the separate instruction manual "NC32 INTRABEAM Floor Stand".

- | | |
|---|---|
| 1 | Axis of motion |
| 2 | Connector panel |
| 3 | Stand base with locking tabs |
| 4 | Control panel |
| 5 | Storage tray |
| 6 | XRS X-ray source
(part of the INTRABEAM Core System) |
| 7 | Handgrip with magnetic brake release button |
| 8 | Head of the NC32 INTRABEAM Floor Stand |

(Blank page, for your notes . . .)

Lined area for notes.

Preparations for use

Requirements for system operation.....	85
Staff requirements	85
Space requirements	85
Electrical requirements.....	86
Visual system check.....	86
Positioning the INTRABEAM cart	87
 Layout of the treatment facility.....	 89
 Connecting	 92
 Workflow	 94
 Powering up the INTRABEAM Core System	 96
Working mode selection.....	98
Application software - main window	100
Menu bar	101
Tab bar.....	105
Data import and export	106
 Equipment Inventory	 112
"XRS X-Ray Sources" group	114
Recalibration of XRS X-ray sources.....	120
"QA Instruments" group	123
"Applicators" group.....	127
 System Quality Assurance.....	 130
Video instruction, video window	136
XRS probe adjustment using PAICH	138
Dynamic Offsets (adjusting the beam deflection)	142
PDA Source Check (isotropy).....	146

PAICH Output Check (dose rate)	153
Only for predecessor systems: checking the ERM external radiation monitor. 161	
Entering treatment parameters	163
Checking the treatment parameters	173

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Staff requirements

**WARNING**

Operation of the INTRABEAM system by untrained staff can damage the system or inadvertently start safety-critical functions and thus generate injury through radiation.

The INTRABEAM system is designed for to be only used by qualified radiotherapists, surgeons, medical physicists and OR staff, who have been properly trained in the use of the INTRABEAM system.

Space requirements

Adequate space to operate the INTRABEAM Core System must be provided both outside the patient care area and for treatment in the treatment room.

The INTRABEAM Core System produces X-rays and its components are sensitive devices.

Vibrations and shocks can cause damage and malfunctions of the sensitive components of the INTRABEAM Core System.

- Make sure that the components of the INTRABEAM Core System are always set up in a vibration-free position.

Strong magnetic fields can affect the XRS X-ray source.

- Do not expose the XRS to strong magnetic fields.

The specified ambient requirements, see Page 259, must be met.

Electrical requirements

The line voltage required and the maximum power consumption to be expected for the complete INTRABEAM system (without INTRABEAM stand) are specified on the type label and in the technical data.

The wall outlet must be provided with a properly connected protective ground conductor.

Do not use power strips; only use single plugs or an isolating transformer.

Visual system check

When starting to prepare the system for use, visually inspect the individual components of the INTRABEAM Core System.

**WARNING:**

If the components of the INTRABEAM system are not checked (visually) prior to use, defective components may be used without being noticed and so result in incorrect measurement. This can lead to incorrect radiation, amongst other things.

- Ensure that the components do not display any obvious damage or changes, such as:
 - deformation, dents, scratches
 - contamination, discoloring
 - kinked, broken or frayed cables.

**WARNING:**

If the user makes use of non-approved components for taking measurements (instead of the IC ionization chamber and UNIDOS E dosimeter), this can lead to incorrect measurements which, in turn, can cause incorrect radiation among other things.

- Check (visually) that only ZEISS-approved components are used.
- Check that the components are calibrated.
- Check that all components are complete and that the right number of components are present.

If scratched or otherwise damaged data carriers (CD or DVD) with calibration data are used, it may not be possible to import the data.

- Prior to use, check that the data carriers are undamaged.

The use of non-approved components for measurement can lead to incorrect measuring results.

- Check (visually) that only ZEISS-approved components are used.

NOTE

If the interlock surfaces are contaminated, the conditions for the interlock may not be met and the INTRABEAM system may not generate X-rays.

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Positioning the INTRABEAM cart

Relocating the unit

- Disconnect the unit from line power and, if necessary, unplug the potential equalization cable from the wall outlet. Properly wrap the power and potential equalization cables around the bolts on the back of the INTRABEAM cart.
- Stow away all single components in the respective transport cases and, if necessary, in the drawer.
- Push in all pull-out shelves and the keyboard tray and close the drawer.
- Release the locking tabs of the casters by pushing up the tabs.
- Push or pull the INTRABEAM cart only on the narrow side using the hand rail provided.
- The casters on the back are casters for straight-ahead travel which can be locked in the movement direction. Operate the tab of one of the casters on the back to lock it in the travel direction. This facilitates straight-ahead movement of the INTRABEAM cart.
- Avoid collisions of any kind and do not move over steps and edges.
- Do not park the INTRABEAM cart on slopes.

**CAUTION**

If the INTRABEAM cart is positioned close to the patient area, germs and infections may be spread.

- During treatment, the INTRABEAM cart must not be positioned close to the patient area, i.e. it must be positioned at least 1.5 m away from the patient.

Securing the INTRABEAM cart in the new position

- Secure the INTRABEAM cart in position by pressing the locking tabs of the front casters. Make sure that the INTRABEAM cart is in a stable, horizontal position and cannot roll away.

Layout of the treatment facility

A suggested layout of equipment and personnel to administer treatment with the INTRABEAM Core System is shown in the figures on Page 90ff.

The treatment room should be regarded as divided into a sterile area and a non-sterile area, as shown in the diagram.

Sterile area

Sterile area refers to an area where a sterile field is maintained, and which includes space for the patient and surgical team, as well as those INTRABEAM Core System components that contact the patient either directly or indirectly.

If non-sterile materials and/or components are brought into the sterile OR area, germs may be spread in this OR area.

- Make sure that only sterile materials and/or components are brought into the sterile OR area.

CAUTION

If no drapes are used, germs may be spread in the sterile OR area.


- Always cover the INTRABEAM stand with a sterile drape when using it in the OR.

An individual in the sterile field is required to handle and operate sterile INTRABEAM Core System components during equipment setup and pretreatment verification procedures. This individual should be familiar with the tasks required to verify proper operation of the INTRABEAM Core System components, and skilled in properly performing those tasks using aseptic techniques.

Layout during tumor bed treatment

e.g. for breast cancer treatment

- 1 OR, sterile area
- 2 OR, non-sterile area
- 3 Patient
- 4 XRS with applicator (sterile)
- 5 INTRABEAM stand (draped)
- 6 Surgical assistant
- 7 Surgeon
- 8 Anesthetist
- 9 Surgical instruments
- 10 OR nurse
- 11 PRS 500 control console
Minimum distance from control console to patient: 1.5 m
- 12 User terminal
- 13 Personnel for radiation therapy (physicist),
operation of control unit and user terminal
- 14 Desk or system cart
In addition to PRS 500 control console (11) and user terminal (12), it accommodates the complete INTRABEAM Core System equipment (stored in the respective transport and storage cases):
 - PDA
 - PAICH
 - INTRABEAM drape
 - Cables
 - Transport and storage cases
 - V-block guide and X-block

Connecting

**WARNING:**

The use of defective connectors or cables can lead to function failure or damage of the INTRABEAM Core System. This may lead to incorrect radiation.

- Visually check all connections and cables before every use of the INTRABEAM Core System.
- If you detect any damage, have it repaired by ZEISS Service.

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
Setting up the INTRABEAM Core System

- Make the following cable connections:
 - Keyboard to user terminal (KB/MOUSE)
 - User terminal to PRS 500 control console, RS232 connection, COM 1 - COM 1
 - User terminal to UNIDOS E dosimeter, RS232 connection, COM 2 - serial
 - PRS 500 control console to XRS X-ray source, color code on PRS 500 control console: green
 - PRS 500 control console to PAICH or PDA, color code on PRS 500 control console: yellow
 - Power supply: connect PRS 500 control console, user terminal and UNIDOS E dosimeter to the power outlet using the respective IEC power plug.
 - Connection for potential equalization
 - External interlock (optional)

Cables - color coding

The INTRABEAM Core System cables are color coded at the end connecting to the PRS 500 control console to match their proper receptacles. They have gray connectors at the end connecting to components in the sterile area.

Color coding		
Cables	Connection to the components	Connection to the PRS 500 control console
XRS	gray	green
PDA/PAICH	gray	yellow

Workflow

If the necessary steps for system verification are performed incorrectly, the XRS probe can be damaged. As a result, the INTRABEAM system can switch off so that radiation can no longer be emitted by the XRS X-ray source.

- Always keep to the specified workflow when working with the INTRABEAM system.


Workflow

Starting the PRS 500 System

First switch on the PRS 500 control console,
then the UNIDOS E dosimeter, if provided,
and then the user terminal.

Check the system components.
For New Component Installation: Data
Entry

System check

Entry of
Treatment Parameters

Checking of
Treatment Parameters

Treatment

Evaluation of
Treatment Parameters

Checking of Treatment Parameters

If a residual dose remains to be adminis-
tered, enter the corresponding parameters
and complete the treatment using the
"Continue with previous patient" function.

Powering up the INTRABEAM Core System

The electrical components of the INTRABEAM Core System, i.e. the PRS 500 control console, UNIDOS E dosimeter and user terminal, are line-operated.

Powering up the electrical components

In the "I" position, the line voltage is switched on.

In the "O" position, the line voltage is switched off.

- Switch on the electrical components of the INTRABEAM Core System in the following sequence:
 - 1 PRS 500 control console
After the PRS 500 has been switched on, it automatically performs a self-test.
 - 2 UNIDOS E dosimeter
After power-on, allow the UNIDOS E dosimeter to warm up for at least 5 minutes before you start the measurements.
 - 3 User terminal

Starting the user terminal

- Switch on the user terminal by pressing the <POWER> control button. This button enables you to switch between the operating mode and standby mode of the user terminal, see Page 58.

Starting the application software


After the user terminal has been started, the start routine of the "PRS 500 Control System" application software is initiated. As soon as the application software has been started, the login window is displayed.

You need a password to be able to use the application software. The login passwords will be handed over to you separately.

Any password changes are the user's own responsibility!

Login window of the application software

- Enter your user ID (user name) (1) and your password (2).

On the initial installation of your INTRABEAM Core System, ZEISS Service will enter the following basic settings:

- Definition of the name and password rules with the user.

User groups

The application software provides five user groups:

Administrator	Physicist	Physician	Radiotherapist	Service
---------------	-----------	-----------	----------------	---------

The physician and physicist user groups must confirm the values checked during verification of the treatment parameters. Make sure that only medical physicists are assigned to the physicist user group and only radiation oncologists to the physician user group.

The administrator can only manage users, but not operate the application software. Every user can create and manage further users with equal rights within his/her user group:

Once the application software has been started, the "Administrative Tools" window for user management can be activated via the following path: Menu bar, <Options> / <Administrative>.

The "Administrative Tools" window opens, and the following functions are available to you in a drop-down menu:

New Account	Overwrite Password	Delete Account
Change Password	List Users	Preferences

For more information, please see Page 103.

Working mode selection


After you have logged in to the application software with your user name and password, you have to select the working mode you wish to use.

- Offline mode, <Work Offline> button (1)
- Online mode, <Start Comm> button (2)

Offline mode, <Work Offline> button (1)

"Work Offline" means working with the user terminal only. This mode can only be used for administrative work. You can:

- manage the components of your INTRABEAM Core System using the "Equipment Inventory" screen mask,
- manage treatment data using the "Review Treatment" screen mask.

Online mode, <Start Comm> button (2)

The Online mode is the full mode and allows you to use all features of the application software. You can:

- manage the components of your INTRABEAM Core System using the "Equipment Inventory" screen mask
- verify your system, i.e. calibrate specific components of your system using the "System Quality Assurance" screen mask
- enter patient data
- enter the treatment parameters for an upcoming treatment session using the "Enter Prescription" screen mask,
- have the treatment parameters checked by the radio oncologist (physician) using the "Verify Treatment Parameters" screen mask,
- conduct and monitor treatment using the "Treatment" screen mask,
- manage treatment data using the "Review Treatment" screen mask.
- The "Calibration" screen mask (Page 210) permits you:
 - synchronize the internal clock of the PRS500 control console with the clock of the user terminal,
 - determine the deflection sensitivity,
 - determine the linearity of the IRM internal radiation monitor,
 - calibrate the air pressure and temperature sensors,
 - only for predecessor systems: adjust the dilator.


Application software - main window

The header of the application software's main window contains the following elements:

- 1 Control buttons for the main window
- 2 Title of the application software
- 3 Menu bar
- 4 Tab bar
- 5 Heartbeat display


Control buttons for the main window

- 1 Minimize the window
- 2 (Disabled)
- 3 Close the main window, i.e. abort the application software

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1 <File>

- <Exit>
This function is used to end the application software.

2 <Reports>

- <Open>

All reports available in the INTRABEAM system are listed.

- Select the required report to have it displayed.

- <Open Latest>



After system verification using the different System Quality Assurance routines, the current results, i.e. the latest SQA reports, can be displayed, printed out or saved as PDF files:

- <PDA> Report of the PDA Source Check
- <PAICH> Report of the PAICH Output Check
- <ERM> (Only for predecessor systems: Report of the ERM Check).

- <Print Latest>

After system verification using the different System Quality Assurance routines, the current results, i.e. the latest SQA reports, can be printed out and saved as PDF files. As soon as you re-run an SQA routine, the data of the previous report file is overwritten.

- <PDA> Report of the PDA Source Check
- <PAICH> Report of the PAICH Output Check
- <ERM> (Only for predecessor systems: Report of the ERM Check).

Creating and saving PDF files


Select the "CutePDF Writer" printer in the printer selection. This "printer" creates a PDF file. The "SAVE" function permits you to save the PDF file to your system, e.g. to D:\intrabeam\sqa-results.

Burning SQA reports to CD or copying them to a USB stick


You can also burn SQA reports created as PDF files to CD, e.g. to print out the reports on a different computer. For more information, please see Page 108.

3 <About>

A window is opened that informs you about the software versions and your INTRABEAM system (serial number).

4 <Tools>


– <Applicator Sterilization>

The menu item <Applicator Sterilization> is used to document how often an applicator has been sterilized.

Several applicators are often provided in the OR for radiation therapy treatment. All these applicators are then considered to be "contaminated" and need to be resterilized.


Since the system is only able to automatically record the applicator selected in "Enter Prescription", once the treatment is started, the sterilization of the other applicators has to be documented manually. We recommend that applicator sterilization be documented right during treatment. For more information, please see "Documenting the number of applicator sterilizations", Page 196.

– <XRS Usage> / <View>


The usage data (operating times) for each X-ray source of your INTRABEAM system is displayed.

The display shows the overall radiation time in hours and the proportion of the treatment time in hours.

– <Administrative>

The "Administrative Tools" window is opened.

The "Administrative Tools" window provides the following options in a pull-down menu (1):

- New Account
You can add further users to your user group.
- Change Password
The user currently logged in can change his/her password.
- Overwrite Password
This function is **only** intended for the administrator.
- List Users
Lists all users logged in to the system.
- Delete Account
Permits a user to be deleted from the user list.
- Preferences
Setting for signature regulation:
If you remove the tick from the entry box, this means that the treatment data no longer needs to be signed by the physician.

If you remove the tick from the entry box, with the result that no second physician needs to check the treatment planning, this may violate existing legal provisions.

– <Set CC volume>

Used for setting the volume of the built-in loudspeaker of the PRS500 control console. The set volume is reset after a restart of the software.

- <Log in>
Reactivates the login window for changing the user.

Tab bar

The tab bar contains all major function modules of the application software arranged from left to right in their logical sequence, i.e. in a regular workflow, you activate these function modules one after the other.

Click on the tab to branch to the associated screen mask.

- Equipment Inventory
Capturing the data of the existing system configuration, see Page 112.
- System Quality Assurance
Verifying the system (SQA), see Page 130.
- Enter Prescription
Entering the treatment parameters, see Page 163.
- Verify Treatment Parameters
Verifying the treatment parameters, see Page 173.
- Treatment
Conducting the treatment, see Page 180.
- Review Treatment
Display of the treatment parameters, see Page 199.
- Calibration
Conducting calibrations, see Page 210.

Heartbeat display

On the extreme right in the tab bar, the heartbeat of the application software is displayed in the form of two boxes changing in color. As long as the heartbeat display moves, the application software is working. The heartbeat display also operates in the offline mode.

Data import and export

Several icons for data import and export are available to you at the computer level of the user terminal.


Data import and export have not been implemented in the application software, but are functions provided at the computer level of the user terminal.

1 Icons for data import

You will find a group of four icons for data import at the top right of the screen.

2 Icons for data export

A group of six icons for data export is provided on the bottom right of the screen.

Importing data

Files can be copied from a CD or DVD to the hard drive of the user terminal.


Vibrations and shocks that affect the user terminal while data is being copied to the hard drive can result in the import of incorrect or incomplete data.

- Make sure that the user terminal is always set up in a vibration-free position.

NOTE

Unintended moving or deletion of data sets can prevent the system verification or treatment from starting.

- Make sure that the INTRABEAM system is used only by radiotherapists, surgeons and OR staff properly trained in the use of the system.

**Importing calibration files of the XRS X-ray source.**

- Insert the CD with the XRS calibration data in the CD drive.
- Click the "XRS_Cal_Files from CD" icon. A directory (window) containing the XRS calibration files on the CD opens.
- Click the "XRS_Cal_Files on Local disk" icon. A window opens.
- Select the required XRS calibration files in the "XRS_Cal_Files from CD" window and copy them to the open "XRS_Cal_Files on Local Disk" window. This way, XRS calibration files can be saved to the hard drive of the user terminal.

Importing calibration files of the applicators:

- Insert the CD with the applicator calibration data in the CD drive.
- Click the "App_Cal_Files from CD" icon. A directory (window) containing the applicator calibration files opens.
- Click the "App_Cal_Files on Local Disk" icon. A window opens.
- Select the required applicator calibration files in the "App_Cal_Files from CD" window and copy them to the open "App_Cal_Files on Local Disk" window. This way, applicator calibration files can be saved to the hard drive of the user terminal.

Saving data**WARNING**

The INTRABEAM system is not intended for permanent data archiving. The calibration and treatment data may be lost or no longer be reproducible if no backup copies are made. This can also lead to over-irradiation.

- Make regular backups.
- All users are responsible for archiving their own data.

Data on the user terminal can be deleted by every subsequent user.

In the event of repairs, all data is deleted in the user terminal. Non-archived data is lost.

You can save data e.g. to a CD or DVD.

If sufficient storage space is no longer available, the INTRABEAM system informs you that files should be archived or deleted, if no longer required.

Exporting data

Files (treatment reports and verification reports) can be copied from the hard drive of the user terminal to CD or a USB stick.

Exporting files to CD

- Insert a blank CD in the CD drive.
- Double-click on the "Copy Files to CD" icon. A window is opened.
- Open the source file window (e.g. report file) by double-clicking on the relevant icon.
- Select the files required and copy them to the "Copy Files to CD" window.

The selected files can be burnt to CD.

- Select the menu item "Write these files to CD" in the Files menu.

Do not close the door of the CD/DVD drive bay during the burning process.

- Close the two windows when the burning process has been completed.

Exporting files to USB stick

Two USB ports are provided in the connector panel on the back of the user terminal.

NOTE

The connection of non-approved devices and the incorrect connection of devices to the user terminal may damage the user terminal and result in the loss of data.

- Do not connect two USB sticks simultaneously.
- Do not connect USB sticks with more storage capacity than 4 GB.
- Do not connect USB hard drives or any other USB devices.

Procedure:

- Insert a USB stick.

- Open the source file window (e.g. report file) by double-clicking on the relevant icon.
- Double-click on the "Copy Files to USB" icon. A window is opened.
- Select the files required and copy them to the USB window.
- Close the two windows when the copying process has been completed.
- Remove the USB stick from the user terminal.

(Icon "Error log on local disk" - currently disabled)

Icon "SQA Results on Local Disk"

A double click on this icon opens a directory (window) containing the SQA results (SQA reports).

The chapter "Verifying the system (SQA)" on Page 134 explains how to export SQA reports from your user terminal.

Icon "Treatment Records"

Saves PDF printouts of the Treatment Records

Support for ZEISS Service

Icon "Inventory on Local Disk"

The "qaini.bin" file is displayed in a separate window. This file provides ZEISS Service with information about the system; it contains all information from "QA Instruments".

Icon "Logfiles on Local Disk"

The "Session*.ibs" file is displayed in a separate window. This file supports ZEISS Service in the event of malfunctions. The files contain information about the relevant session.

The "qaini.bin" and "Session*.ibs" files provide ZEISS Service with information about the system. This information permits ZEISS Service to perform malfunction analysis, to undertake remote diagnosis, if possible, or to prepare service calls.

If required, ZEISS Service will ask you to send these two files to them by email.

(Blank page, for your notes . . .)

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. A vertical pink or light red margin line runs down the right side of the page. The paper appears to be from a notebook or a standard composition book. There are no markings, text, or drawings on the page.

Equipment Inventory

"Equipment Inventory" screen mask

The "Equipment Inventory" screen mask allows you to capture, display and update the data of the individual components of your system. You can also monitor the periodic calibrations required.

The relevant histories give you an overview of the calibrations and the persons conducting them.

Activating the function

- Activate the "Equipment Inventory" screen mask by clicking the <Equipment Inventory> tab.

The following component groups are listed in the "Equipment Inventory" screen mask:

- XRS X-ray sources
- Quality assurance instruments
- Applicators

To manually edit the data of a group, press the associated <Edit> button.

"Equipment Inventory" screen mask

- 1 "XRS X-Ray Sources" group
Use this table to manage all XRS X-ray sources of your INTRABEAM Core System, for details see Page 114.
- 2 "QA Instruments" group
Use this table to manage all verification devices of your INTRABEAM Core System, for details see Page 123.
- 3 "Applicators" group
Use this table to manage all applicators of your INTRABEAM Core System. For details see Page 127.

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | Review Treatment | Calibration

1

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | Review Treatment | Calibration

X-Ray Sources Edit

XRS Name	XRS No.	Cal kV	Cal by	Cal Date	Cal Due	Verified by	Status	O.C. (Gy/min)	Check Date	Interval	Notes
1st XRS 40	000170	40.0	Nrm	04-Apr-07	Annual	bess	Active	1.877	13-Jun-07	Every Treatment	1st XRS 40
1st XRS 50	000170	50.0	zobess	04-Apr-07	Annual	bess	Active	2.849	12-Jun-07	Every Treatment	1st XRS 50

2

QA Instruments Edit

Instrument	Model	Serial No.	Cal Factor	Cal Units	Cal Date	Cal Due	Calibrated by	Verified by	Status	Notes
Control Console	CC500	0001	N/A		01-Jan-07	Permanent	CZ	Bess	Active	1st CC
Dosimeter	UNIDOS-E	00250	1.000E+0		21-Apr-07	Annual	PTW	Bess	Active	1st Dosim.
ERM Test Adapter	ETA400	00002	N/A		01-Jan-07	Permanent	CZ	Bess	Active	1st ETA
ERM	ERM400	00017	N/A		01-Aug-06	Permanent	CZ	Bess	Active	1st ERM
Ion Chamber	PTW23342	N23342-1006	1.141E+9	Gy/C	13-Sep-06	Annual	CZ	Bess	Active	1st IC
PAICH	PAICH450	507190	-2.914E+1	Deg C	01-Jan-07	Permanent	CZ	Bess	Active	1st PAICH
PDA	PDA400	00038	N/A		01-Mar-06	Permanent	CZ	Bess	Active	1st PDA

3

Applicators Edit

App. Name	Type	Size	Serial No.	Cal by	Cal Date	Verified by	No. of Sterilizations	Status	Notes
4.0 Appl.	Spherical	4.0	264704	zobess	02-Mar-07	Physicist	0 of 100	Active	1st 4.0 Appl.
4.5 Appl.	Spherical	4.5	264804	zobess	02-Mar-07	Physicist	0 of 100	Active	1st 4.5 Appl.
5.0 Appl.	Spherical	5.0	264904	zobess	02-Mar-07	Physicist	0 of 100	Active	1st 5.0 Appl.
1.5 Appl.	Spherical	1.5	265101	zobess	02-Mar-07	Physicist	0 of 100	Active	1st 1.5 Appl.
2.0 Appl.	Spherical	2.0	265201	zobess	02-Mar-07	Physicist	0 of 100	Active	1st 2.0 Appl.
2.5 Appl.	Spherical	2.5	265301	zobess	02-Mar-07	Physicist	0 of 100	Active	1st 2.5 Appl.
3.0 Appl.	Spherical	3.0	265401	zobess	02-Mar-07	Physicist	0 of 100	Active	1st 3.0 Appl.
3.5 Appl.	Spherical	3.5	265501	zobess	02-Mar-07	Physicist	0 of 100	Active	1st 3.5 Appl.

"XRS X-Ray Sources" group

Use this group to manage all X-ray sources of your INTRABEAM Core System.

The XRS X-ray radiation source data contained in your INTRABEAM Core System delivery package is the ZEISS calibration data. This data is installed by ZEISS Service staff during the initial installation of your INTRABEAM Core System.

You can enter further data manually, see "Data import and export" Page 106.

- To edit the data of this group manually, press the relevant "Edit" button.

An editing window opens.

Control functions

Editing individual boxes:

- Directly click on the box you want to edit.

This may be a menu box or an input box.

When you click on a menu box, a pull-down menu is opened which offers you all setting options for selection or shows all existing entries (e.g. numerical values).

When you click on an input box, you can directly enter numerical values.

Scroll function:

Some of the boxes have a button icon on their left (1). You can use this button icon to scroll forward and backward in the setting options provided in the relevant box.

Button functions:

- 1 <Enter> button
Saving the current changes. You can go to another line after <Enter>.
- 2 <Continue> button
Saving all current changes and returning to the "Equipment Inventory" screen mask.
- 3 <Cancel> button
Rejecting all current changes and returning to the "Equipment Inventory" screen mask.

Inactive

✓ Active

1

"XRS Inventory" screen mask

[illegible]

Table columns

- 4 XRS name
Name of the XRS X-ray radiation source, any name can be entered
- 5 XRS nr.
Serial number of the XRS X-ray radiation source, contained in the calibration file and indicated on the XRS X-ray radiation source
- 6 Cal kV
50 kV, 40 kV: The voltage determines the maximum X-radiation output.
- 7 Calibrated by
The name is entered automatically
- 8 The calibration date is entered automatically
- 9 Calibration due (calibration interval)
Observe the legal requirements for the XRS X-ray radiation source; we recommend calibration every year.

The calibration interval of the XRS X-ray radiation source commences with the commissioning of your INTRABEAM system. As there can be a period of several months between the calibration and the commissioning of the system, "Biennial" (calibration every two years) has been entered in this field at the time of the initial delivery. You will however automatically be made aware of the expiry of the calibration interval by ZEISS within one year of commissioning.

 - Reset the calibration interval to "annual" after expiry of the first calibration interval .
- 10 Verified by
Name of the person loading and accepting the calibration file.
- 11 Status
"Active" or "Inactive"
- 12 "DR" (dose rate), "Check date" and "IC Check Interval"
The OC" (Gy/min), "Check Date" and "Interval" boxes contain entries which are the result of the "PAICH XRS Check" of the system
- 13 IC test interval
Test interval of the PAICH XRS check (dose rate) using the ionization chamber IC
- 14 Notes; individual remarks possible

"XRS Inventory" screen mask

XRS Name	Status	Calibration Due	IC Check	Interval	Notes
1st XRS 50	Active	Annual	Every T	reatment	1st XRS 50

Buttons:

Calibration intervals "Cal. interval"

The "Cal Due" column indicates the calibration intervals of the relevant XRS X-ray source.

You have to enter the calibration intervals for each XRS manually in the "Cal Due" editing box.

ZEISS recommend that you have the XRS X-ray radiation source calibrated annually .

- Please observe the calibration regulations on Page 209.
On expiry of the calibration interval, the XRS X-ray radiation source will automatically be set to "Inactive".

Status of an X-ray source

The "Status" column indicates the status of the relevant XRS X-ray radiation source.

Active: The component is available for verification.

Inactive: The component is not available for verification.

The "Status" editing box permits you to change the status of an XRS X-ray radiation source from "Inactive" to "Active".

If you set an XRS X-ray radiation source whose calibration has expired to "Active", you are requested to confirm this with your password, i.e. continued operation is at your own responsibility.

This setting to "Active" must be repeated after every software restart.

(Blank page, for your notes . . .)

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. A vertical red margin line runs down the right side of the page. The paper appears to be from a notebook or a standard composition book. There are no markings, text, or drawings on the page.

Recalibration of XRS X-ray sources

We recommend that you have the XRS X-ray radiation source calibrated annually.

Your XRS X-ray radiation source is calibrated in the ZEISS manufacturing plant. In the calibration process, the characteristic of the XRS X-ray radiation source, the dose rate [Gy/min] relative to the treatment depth [mm] is measured at 40 kV and 50 kV. This is done in a water phantom, a special container filled with water.

The calibrated XRS X-ray radiation sources will be returned to you with a calibration file. You must copy this calibration file to the application software on your user terminal, see the chapter "Data import and export" Page 106.

- Copy this calibration file to the software of your INTRABEAM Core System in the following folder:
D: \ intrabeam \ cal \ xrs.
- The old calibration file is automatically transferred to the archive directory.

Saving the calibration file


The application software automatically recognizes new calibration data. The next time the software starts, the application software will ask you: "Please verify the XRS calibration data".

If you acknowledge with OK, a window with the new calibration file is opened.

You will see the same screen mask later when you select an XRS X-ray radiation source (line in the table) with a double-click. In this case, the screen mask is displayed for viewing only.

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | Review Treatment | Calibration | 

This image shows a single page of blank, lined notebook paper. The paper is white and features evenly spaced, horizontal blue lines running across its entire width. There are no vertical margin lines or other markings present. The lighting is uniform, and the lines are clearly visible against the white background.

"QA Instruments" group

Use this group to manage all verification components of your INTRABEAM Core System:

- PRS 500 control unit - UNIDOS E Dosimeter
- Ionization chamber (IC) - Photodiode array (PDA)
- PAICH probe adjuster/ionization chamber holder

The data of the components contained in your INTRABEAM Core System delivery package is entered in this group by ZEISS Service during the initial installation of your INTRABEAM Core System. Further data can be manually entered by you.

- See the screen mask "Component management" Page 112: To edit the data of this group manually, press the relevant "Edit" button. The relevant editing window is opened.

WARNING:

If the operator does not observe the calibration intervals specified for the IC ionization chamber, the measurements performed during the system check can be incorrect. This can lead to over or under-irradiation.

- Have your IC ionization chamber calibrated every two years in accordance with the recommendations of PTW.

Control functions

For a description of how to edit the individual entries, see Page 114.

Button functions

- <New> Adding a new component to the table, i.e. a new line is created in the table.
- <Enter> Saving the current changes. You can go to another line after <Enter>.
- <Delete> Deleting the item currently selected, i.e. the line currently selected in the table is deleted.

Table columns and control functions

For a description of how to edit the individual entries, see Page 114.

- 1 "Instrument" column (devices of your INTRABEAM Core System)
If your INTRABEAM Core System includes several devices of the same type (e.g. dosimeter), please note that only one of these devices can be "Active" at a time.
- 2 "Serial no." column (serial number of the relevant device)
The system requires the serial number of the device for communication and reports.
- 3 "Cal Factor" column (calibration factor)
The calibration factor is provided to you by the organization performing the calibration in the calibration file for the relevant device.
Calibration factor for IC ionization chamber:
The calibration factor with the measuring unit Gy/C is used for dose rate measurement in the "PAICH Output Check " screen mask, see Page 157.
PAICH data: the data of the PAICH probe adjuster/ionization chamber holder is taken from the result of the temperature sensor calibration procedure, see Page 223. This data is automatically entered here after completion of the calibration.
- 4 "Cal Due" column (calibration interval)
Please note the applicable regulations for the different devices.
- 5 "Status" column
Active: The device is available for further procedures.
Inactive: The device is no longer available.
- 6 <Delete> button
Deletes the item currently selected, i.e. the line currently selected in the table is deleted.
- 7 <New> button
Adds a new component to the table, i.e. a new line is created in the table.
- 8 <Enter> button
Saves the current changes. You can go to another line after pressing <Enter>.

[illegible]

(Blank page, for your notes . . .)

Blank lined paper for writing.

"Applicators" group

Use this table to manage all applicators of your INTRABEAM Core System.

Each applicator is supplied to you with a calibration file, see the chapter "Data Import and export" Page 106.

- Copy these files to the software of your INTRABEAM Core System using the following folder:
D:\Intrabeam\cal\applicators.

Control functions

Editing individual boxes:

- Directly click on the box you want to edit.

This may be a menu box or an input box.

When you click on a menu box, a pull-down menu is opened which offers you all setting options for selection or shows all existing entries (e.g. numerical values).

When you click on an input box, you can directly enter numerical values.

Scroll function:

Some of the boxes have a button icon on their left (1). You can use this button icon to scroll forward and backward in the setting options provided in the relevant box.

Button functions:

- | | |
|--|---|
| <div data-bbox="295 1079 456 1126" data-label="Text">Inactive</div> <div data-bbox="295 1137 419 1178" data-label="Text">Active</div> <div data-bbox="242 1279 258 1303" data-label="Text">1</div> | <p><Enter> Saving the current changes. You can go to another line after <Enter>.</p> <p><Continue> Saving all current changes and returning to the "Equipment Inventory" screen mask.</p> <p><Cancel> Rejecting all current changes and returning to the "Equipment Inventory" screen mask.</p> |
|--|---|

"Applicators Inventory" screen mask

- 1 Name (editing box)
This is an editing box for any type of entry; you can enter a personal text (name) here.
- 2 Status (editing box)
This editing box permits you to change the status of the relevant applicator from "Inactive" to "Active".
Active: The applicator is available for further treatment.
Inactive: The applicator is not available for further treatment.

If you set an applicator that has reached the maximum admissible number of sterilizations back to "Active", you are requested to confirm this with your password, i.e. continued operation is at your own responsibility.

**WARNING**

Exceeding the maximum admissible number of sterilizations (including cleaning and disinfection) can damage the applicator and lead to incorrect radiation.

- Make sure that all sterilization cycles are counted and that no spherical applicator is used that has been sterilized more than 100 times.

- 3 Notes (editing box)
This is an editing box for any type of entry; you can enter a personal text (notes) here.
- 4 No. of sterilizations (display box)
The number of sterilization processes is incremented automatically and displayed here. When the specified maximum number of applications has been reached, the applicator is automatically set to "Inactive". A spherical applicator may be sterilized only a limited number of times, and must no longer be used after that.
Please also note the chapter "Documenting the number of applicator sterilizations", Page 196.
- 5 Status (display box)
The "Active" or "Inactive" status of the relevant applicator is indicated.

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | Review Treatment | Calibration

1	2	3	4	5
Name	Status	Notes		
4.0 Appl.	Active	1st 4.0 Appl.		

App. Name	Type	Size	Serial No.	Cal by	Cal Date	Verified by	No. of sterilizations	Status	Notes
4.0 Appl.	Spherical	4.0	264704	zobess	02-Mar-07	Physicist	0 of 100	Active	1st 4.0 Appl.
4.5 Appl.	Spherical	4.5	264804	zobess	02-Mar-07	Physicist	0 of 100	Active	1st 4.5 Appl.
5.0 Appl.	Spherical	5.0	264904	zobess	02-Mar-07	Physicist	0 of 100	Active	1st 5.0 Appl.
1.5 Appl.	Spherical	1.5	265101	zobess	02-Mar-07	Physicist	0 of 100	Active	1st 1.5 Appl.
2.0 Appl.	Spherical	2.0	265201	zobess	02-Mar-07	Physicist	0 of 100	Active	1st 2.0 Appl.
2.5 Appl.	Spherical	2.5	265301	zobess	02-Mar-07	Physicist	0 of 100	Active	1st 2.5 Appl.
3.0 Appl.	Spherical	3.0	265401	zobess	02-Mar-07	Physicist	0 of 100	Active	1st 3.0 Appl.
3.5 Appl.	Spherical	3.5	265501	zobess	02-Mar-07	Physicist	0 of 100	Active	1st 3.5 Appl.

Cancel Continue

System Quality Assurance

System Quality Assurance (SQA) is the verification of the INTRABEAM system.

If System Quality Assurance has been successfully completed, the green indicator lamp of "SQA Valid" is lit. If the steps for system verification are performed incorrectly, the INTRABEAM system no longer permits radiation to be emitted by the respective XRS X-ray source.

- For System Quality Assurance, at least the PDA Source Check and the PA-ICH Output Check must be successfully completed.

The described verification procedures must be performed prior to actual irradiation of a targeted lesion with the INTRABEAM Core System.

NOTE

If the interlock surfaces are contaminated, the interlock may not be satisfied and the INTRABEAM system may not generate X-radiation.

- Check (visually) that the interlock surfaces are clean prior to every use.

Validity of verification

Verification was successful if the green indicator lamp is lit. The verification of the system is now valid for 36 hours maximum. Time is now counted backward, starting at 36 hours.

The time left for which the verification still remains valid is displayed in the window next to the indicator lamp. The time display is updated every time an XRS X-ray source and the associated voltage are selected.


The verification ceases to be valid after 36 hours or when the treatment has been completed. The window next to the indicator lamp then displays "0.0" (h).

If the PDA test is successfully completed again within 36 hours after successful system verification, the PAICH test is not automatically set to "Mandatory" and the green indicator lamp remains lit.

"System Quality Assurance" screen mask

The "System Quality Assurance" screen mask allows you to verify your system by calibrating specific system components.

- Open the "System Quality Assurance" screen mask by clicking on the "System Quality Assurance" tab.

"Choose XRS"

Here you can select the XRS for the subsequent treatment session using its serial number and voltage (kV). The selection box only contains the XRS X-ray sources that have been set to "Active" in the "Equipment Inventory" (Page 114).

"Mandatory"

System verification is obligatory for functions where "Mandatory" has been selected in the system (ticked checkbox).

"Recommended"

For functions where "Recommended" has been selected (ticked checkbox), calibration is not obligatory but can optionally be performed.

Equipment Inventory | System Quality Assurance

Choose XRS: 000170 40.0 kV | IRM Rate [Hz]: 0

Interlock Status: SQA Valid: 0.0 [h]

PAICH/ETA


Probe Adjuster	
Dynamic Offsets	
PDA Source Check	
PAICH Output Check	
ERM Check	

Color code

The color code refers to each individual System Quality Assurance routine. The status of the selected XRS X-ray source (serial no.) for a specific high voltage [kV] (see the "Choose XRS" box in the SQA screen masks) is indicated by colored squares next to the individual verification functions.

- Yellow** The selected XRS X-ray source is available in the system for calibration for this procedure.
- Green** The selected XRS X-ray source has been successfully calibrated for this procedure.
- Red** Error, the selected XRS X-ray source cannot be used until it has passed this verification step, see "Workflow" on Page 94.
- SQA active** In the SQA screen masks, the color code indicates the status of a specific XRS X-ray source (serial no.) for a specific high voltage [kV].
If the color code square is green, this means that the relevant XRS X-ray source is available for use at a specific high voltage [kV] for treatment planning (XRS@kV). For this, at least the color code squares of the PDA and PAICH procedures must be green to indicate that calibration has been successfully completed.

Opening and printing current SQA reports

After system verification using the different SQA routines, you can use the "Reports" menu to display or print out the current results, i.e. the latest SQA reports, for the following verifications:

- PDA PDA Source Check
- PAICH PAICH Output Check
- ERM (only for predecessor systems: "ERM Check")

Selecting an SQA report. <Open>

All SQA reports available in your INTRABEAM system are listed.

- Click on <Open> in the <Reports> menu.
- Select the required SQA report.

Selecting the current SQA report. <Open Latest>

- Click on <Open Latest> in the <Reports> menu and select the SQA routine required: <PDA>, <PAICH> or <ERM>. The relevant latest SQA report is displayed as a table in a separate window.

Use the <Print> button to print out the table and export it as a PDF file.
Use the <Done> button to close the window again.

Printing and saving SQA reports as PDF files

If you wish to export an SQA report for display and/or printout on a different computer, you should print out and save the relevant SQA report as a PDF file.

Recommended storage location: D:\SQResults

To enable clearer subsequent assignment of the PDF file to the treatment session and to facilitate identification among several different SQA reports, we recommend a "text_date_time.pdf" structure for the name of the PDF file. For example:

- PDA_MM_DD_YYYY_hh_mm.pdf

How to create and save PDF files is described on Page 102.

Exporting SQA results

The SQA results can be exported to a USB stick or CD. Export to a USB stick is the faster and easier method.

- At the computer level of the user terminal, double-click on the "SQA Results on local disk" icon (see Page 108). A directory (window) containing the SQA results is opened.

Data export to USB stick

- Insert a USB stick into one of the USB ports on the back of the user terminal.
- Double-click on the "Copy Files to USB" icon. A window is opened.
- Select the SQA results required and copy them to this window.

Data export to CD

- Insert a blank CD in the CD drive.
- Double-click on the "Copy_Files to CD" icon. A window is opened.
- Select the SQA results required and copy them to this window.
- The SQA results required can be burnt to CD. In the "Files" menu, select the option "Write these files to CD".

(Blank page, for your notes . . .)

[illegible]

Video instruction, video window

When you have started an SQA or system calibration procedure, the steps to be performed are shown in a video instruction with displayed text. The application software checks internally whether all verification devices required have been connected and are in position. If this is not the case, the SQA or system calibration procedure cannot be continued.

Comments on the video window

The example shown in the illustration is the video window for the SQA function "PDA Source Check".

The field on the right (1) displays the description (text) for this procedure. The video sequence is shown in the black field on the left (2).

On the bottom left, the "Optical Interlock Status" (3) is displayed (here: PDA photodiode array).

Next to it, you find the status display for connected cables and interlock statuses to be fulfilled, e.g.:

- XRS cable connected
- PDA cable connected
- External Interlock Status

If the associated status indicator (4) lights up green, the specified requirement is met.


- Click on <Continue> (5) to proceed.

If the status indicator is not lit or if the interlock status has not been fulfilled, <Continue> button (5) is disabled. Use <Cancel> button (6) to abort the selected SQA or calibration procedure.

Overview of the video instruction screen mask

- 1 Text display window
This field contains instructions for use in the form of text concerning the currently selected screen mask.
- 2 Video display window
This field contains instructions for use in the form of a video sequence concerning the currently selected screen mask.

- 3 Device indicator
Status display of the optical interlock (in this example: PDA photodiode array).
- 4 Status indicators
If the associated status indicator lights up green, the specified requirement is met.
- 5 <Continue> button
You can only continue with the procedure if the status indicators are lit and the interlock status has been satisfied.
- 6 <Cancel> button
This button permits you to abort the selected SQA or calibration procedure.

XRS probe adjustment using PAICH

If the PDA Source Check was successful, the XRS probe adjustment procedure using PAICH is optional.

You need to perform the XRS probe adjustment if the isotropy verification (PDA Source Check, Page 146) was unsuccessful or the XRS probe is suspected to be bent.

If the XRS X-ray source has already been successfully verified, and if the XRS probe has been subsequently centered in the probe adjustment procedure, it is absolutely vital that you repeat the "Dynamic Offset", "PDA Source Check" and "PAICH Output Check" tests.

Before starting the XRS probe adjustment procedure, it is necessary to do the following:

If you have a PAICH with a micrometer screw:

- Turn out the micrometer screw until the "3" marking on the micrometer scale is visible. This ensures that the micrometer will not damage the XRS probe.
- Insert the XRS probe into the PAICH using the V-block guide.

WARNING

The XRS probe can be mechanically damaged or become unusable when the XRS X-ray source is inserted into the verification devices. This may also lead to incorrect irradiation.

- Always use the V-block guide when inserting the XRS into the verification devices.



WARNING

The XRS probe can be mechanically damaged or become unusable when the orientation of the V-block is incorrect during insertion of the XRS into the verification devices. This may also lead to incorrect irradiation.

- One side of the V-block is designed specifically to hold the XRS X-ray source, while the other side holds the verification devices. This is marked by an appropriate engraving on the V-block.

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | Review Treatment | Calibration


During setup, the XRS can topple over and be damaged or injure the user.

- Connect the PAICH and the XRS to the control console using the appropriate cables.
- Place the XRS on end in the X-block, with the XRS vertical. This position allows the PAICH to be rotated.

Adjusting the XRS probe

Activating the function

- In the "System Quality Assurance" screen mask, select the option "Probe Adjuster".

Video instruction

Press the  button to start an SQA or system calibration procedure, and the application software will visualize the steps to be performed in a video instruction with displayed text, see Page 136. The application software checks internally whether all verification devices required have been connected and are in position. If this is not the case, the procedure cannot be continued.

After completing the installation of the measuring setup, allow it to adjust to the ambient temperature for at least 1 minute before starting measurement.

Measurement and setting of the straightness of the XRS probe is incomplete if the PAICH is not rotated by a full 360°.

- When measuring, rotate the PAICH completely i.e. $\geq 360^\circ$ around the XRS probe being measured.
- Press the <Zero> button (3) in the "PAICH: Probe Adjuster" screen mask, and slowly turn the PAICH through at least 360°.
Hold the PAICH at its base. Do not change your grasp position too frequently when turning.
- The precise total runout is now given in the "Runout" window.
- Rotate the PAICH until the deviation bars are at their highest value.
- Carefully depress the plunger to straighten the XRS probe, then release the plunger. Watch the deflection of deviation bar (2).
- Press the <Zero> button (3) to reset the bars to 0 mm.
- Repeat this procedure until the total runout value (display field "Total") is less than 0.10 mm, recommendation: approx. 0.07 mm.
- Press the <Finished> button (4) to accept the result of XRS probe straightening and exit the routine.

If the runout remains greater than 0.10 mm, all the steps must be repeated. Further equipment verification steps ("Dynamic Offset" Page 142) cannot be performed without successful completion of the XRS probe adjustment procedure.


If you do not succeed in centering the XRS probe to a runout less than 0.10 mm this verification step has not been passed. Please contact ZEISS Service.

The XRS probe can be mechanically damaged when the XRS X-ray source is removed from the verification device.

- Always use the V-block guide when removing the XRS X-ray source from the verification instrument.

Screen mask overview

- 1 Selecting the verification routine, "Probe Adjuster" menu option
- 2 Deviation bars
Graphic display of the XRS probe's deviation from straightness.
- 3 Zero button (reset)
Resetting the deviation bars to zero
- 4 <Done> button (accept)
Accepting the test results or canceling probe adjustment
- 5 User prompt, text display of the current step in the procedure
- 6 Graphic display of the XRS probe concentricity
No function in this screen mask

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When performing the beam deflection adjustment, you align the direction of the electron beam with the mechanical center of the XRS probe. If beam deflection is not adjusted, the electron beam may not fully hit the target after realignment of the XRS probe (probe adjustment) and the verification steps may therefore not be completed successfully. The mechanical alignment of the XRS probe, i.e. centering of the XRS probe (Probe Adjuster function), must always be followed by the electronic alignment of the XRS probe (Dynamic Offsets function).

- Perform this verification step after centering of the XRS probe (Probe Adjuster) or if you imported a new calibration file as a result of recalibration. Otherwise, this verification step is optional.
- Insert the probe into the PDA photodiode array using the V-block guide.

WARNING

The XRS probe can be mechanically damaged or become unusable when the XRS X-ray source is inserted into the verification devices. This may also lead to incorrect irradiation.

- Always use the V-block guide when inserting the XRS into the verification devices.



WARNING

The XRS probe can be mechanically damaged or become unusable when the orientation of the V-block is incorrect during insertion of the XRS into the verification devices. This may also lead to incorrect irradiation.

- One side of the V-block is designed specifically to hold the XRS X-ray source, while the other side holds the verification devices. This is marked by an appropriate engraving on the V-block.

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | Review Treatment | Calibration

- To prevent the XRS from toppling over, place it on end in the X-block, with the XRS vertical.
- Connect the PDA and the XRS to the control console using the appropriate cables.
- Align the +X and +Y markings on both the XRS and the PDA.

Activating the function


- In the "System Quality Assurance" screen mask, select the option "Dynamic Offsets".

Video instruction

Press the  button to start an SQA or system calibration procedure, and the application software will visualize the steps to be performed in a video instruction with displayed text, see Page 136. The application software checks internally whether all verification devices required have been connected and are in position. If this is not the case, the procedure cannot be continued.

- Press the <Start> button. The adjustment is performed automatically for the selected XRS X-ray source at a specific high voltage (XRS@kV). The procedure will take a few minutes.

For other accelerating voltages, the procedure must be performed separately.

When the beam deflection adjustment has been successfully completed, the "Dynamic Offset" color code display changes to green, i.e. you can perform further verification step (PDA X-ray source check Page 146).

If the beam deflection cannot be successfully adjusted, an error message informs you of the reason. Inform ZEISS Service in this case. If the PDA and PAICH color code display already indicated green, you must repeat the PDA and PAICH procedures.


The XRS probe can be mechanically damaged when the XRS is removed from the verification device.

- Always use the V-block guide when removing the XRS from the verification device.

Screen mask overview

- 1 Selecting the verification step
- 2 Display of offset values X, Y
- 3 Stopping the verification step
- 4 Text display of the current step in the procedure
- 5 Starting the verification step
- 6 Graphic display of the beam deviation from the center

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | Review Treatment | Calibration

PDA Source Check (isotropy)

This verification step must always be performed.

The PDA Source Check function is used to verify and, if necessary, adjust the XRS probe's X-ray emission pattern to make it as spherical as possible.

The Isotropy Adjust procedure is performed using maximum beam current and exactly the accelerating voltage selected for the relevant XRS X-ray source.

The system requires that a treatment be performed only at an accelerating voltage setting at which the PDA Source Check test was successfully completed. Performing the PDA Source Check procedure invalidates the results of the subsequent verification test (PAICH), requiring this to be repeated before a treatment can be performed.

It is advantageous to perform the PDA Source Check procedure for an XRS X-ray source multiple times at whatever accelerating voltages are anticipated for treatment before proceeding to the remaining verification tests. This avoids the need to repeat the tests should the treatment accelerating voltage ultimately differ from the one(s) used during isotropy adjustment.

Isotropy adjustment - PDA Source Check

You have to use the PDA for isotropy adjustment. For operation and display, use the "System Quality Assurance/PDA Source Check" menu.

The procedure to adjust the X-ray field of the XRS probe has two major steps:

- Setting of the following operating parameters:
Accelerating voltage, beam current and run time
- Automatic optimization of the X-ray field isotropy.
- Before starting the adjustment procedure, insert the XRS probe into the PDA using the V-block:

[Equipment Inventory](#)[System Quality Assurance](#)[Enter Prescription](#)[Verify Treatment Parameters](#)[Treatment](#)[Review Treatment](#)[Calibration](#)**WARNING**



The XRS probe can be mechanically damaged or become unusable when the XRS X-ray source is inserted into the verification devices. This may also lead to incorrect irradiation.

- Always use the V-block guide when inserting the XRS into the verification devices.

**WARNING**

The XRS probe can be mechanically damaged or become unusable when the orientation of the V-block is incorrect during insertion of the XRS into the verification devices. This may also lead to incorrect irradiation.


- One side of the V-block is designed specifically to hold the XRS X-ray source, while the other side holds the verification devices. This is marked by an appropriate engraving on the V-block.

- To prevent the XRS from toppling over, place it on end in the X-block, with the XRS vertical.
- Connect the PDA and the XRS to the control console using the appropriate cables.

If the user does not correctly position the PDA on the XRS probe, isotropy cannot be adjusted.

- Align the +X and +Y markings on both the XRS and the PDA.

Activating the function

- In the "System Quality Assurance" screen mask, select the option "PDA Source Check".

Video instruction

Press the  button to start an SQA or system calibration procedure, and the application software will visualize the steps to be performed in a video instruction with displayed text, see Page 136. The application software checks internally whether all verification devices required have been connected and are in position. If this is not the case, the procedure cannot be continued.

- To start the isotropy adjustment procedure, press the <Start> button in the "System Quality Assurance / PDA Source Check" menu.

Screen mask for isotropy adjustment

Once you have pressed the <Start> button, X-rays are produced and the software of the PRS 500 control console runs an automatic isotropy optimization process for the XRS probe which will take a few minutes.


The "Photodiode Array [V]" (2) bar graph is displayed in the "PDA Source Check" screen mask, with bars corresponding to the five photodiodes in the PDA.

These bars are labeled +X, -X, +Y, -Y and +Z, and indicate the voltages currently present at the photodiodes.

After the procedure, the four bars +X, -X, +Y, -Y indicate the voltages for optimum isotropy adjustment and the tolerance bars. Because the tolerance bars are displayed relative to the +Z bar, the height of the +Z bar remains fixed at a deviation of 0 V. When the X-ray field is isotropic, the percent deviation of each of the first four bars will be close to zero. The tolerance bars on each of the bars +X, -X, +Y, -Y indicate whether the isotropy is within tolerance or not.

Screen mask overview

- 1 "PDA Source Check" menu item
Selecting the verification step
- 2 Bar graph of sensor signals
- 3 Stopping the procedure
- 4 Text display of the current step in the procedure
- 5 Graphic display of the IRM count rate (counts/time)
- 6 Starting the procedure
- 7 Graphic display of the XRS probe concentricity

Isotropy adjustment out of tolerance

If the tolerance limits have been exceeded, the color code box is red. You have to recenter the XRS probe. Change to the "Probe Adjuster" function for this purpose, Page 138.

Isotropy adjustment within tolerance

If the Automated Isotropy Adjust procedure has been successful, the height of all five bars is within the white tolerance bars, and the color code box is green.

The subsequent "PDA Source Check" window, i.e. the "IRM Rate History" window, opens automatically (see Page 151). The inquiry "Update History and Calibration Files?" is displayed which you can answer using the <Yes> or <No> button.

If you answer with <Yes>, the isotropy data and the IRM rate of the XRS probe are saved and the isotropy adjustment procedure is closed. The IRM rate is subsequently used for treatment planning. You can now perform the "PAICH Output Check", see Page 153. If you answer with <No>, the results are discarded.

The "Reports" menu permits you to print out the isotropy data, see Page 152.

Check the XY diagram on the bottom left of the screen mask which indicates the electron beam runout in the X and Y directions. The red circle represents the absolute limit value, i.e. the system is shut down when this value is exceeded.

In practical use, no more than 50% of this tolerance should be utilized. Therefore, you should recenter the XRS probe when approx. 50% of the admissible maximum runout are reached. Activate the "Probe Adjuster" function for this purpose, see Page 138.

The XRS probe can be mechanically damaged when the XRS is removed from the verification device.


- Always use the V-block guide when removing the XRS from the verification device.

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | Review Treatment | Calibration | 


"IRM Rate History" window, subsequent window to "PDA Source Check"

The diagram shows the time variance of the IRM count rate (counts/time) to the first value and the percentage deviation from the first value.

- 1 <Yes> button
Values are saved.
- 2 <No> button
Values are discarded.
- 3 Display of percent deviation
Percent deviation of the current measurement from the first value ((Day 0)

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | Review Treatment | Calibration

Current report for the "PDA Source Check"

After system verification using the "PDA Source Check" routine, you can display or print out the current result, i.e. the latest SQA report of the "PDA Source Check" routine:


Displaying reports, <Open Latest>

In the <Reports> menu, click on <Open Latest> and select <PDA>. The report of the latest "PDA Source Check" is displayed as a table in a separate window.

Exporting reports of the "PDA Source Check", see Page 108.

The <Print> button permits you to print out the table as a PDF file and export it.

Use the <Done> button to close the window.

PAICH Output Check (dose rate)

This verification step must always be performed.

The "PAICH Output Check" procedure is used to determine the dose rate of the XRS X-ray source by means of the IC ionization chamber and the UNIDOS E dosimeter. The dose rate is checked against a specific, expected value from the calibration data. A correction factor is calculated for the difference between the measured dose rate and the data from the calibration file and is taken into account in subsequent treatment planning.

The isotropy adjustment procedure (PDA Source Check) must have been successfully completed. The successful PDA Source Check is indicated by the green color code of the menu item and the enabled "PAICH Output Check" procedure.

Test setup

This procedure is used to check the dose rate currently being generated by the XRS probe against an expected value in the Z-axis. In order to perform this operation, the PAICH must be attached to the XRS, the IC ionization chamber inserted into the holder and connected to the dosimeter. The application software has already set the UNIDOS E dosimeter to the "Current", "Low pA" measuring mode.

If you wish to change the dosimeter setting at this time, you have to end the application software and make the UNIDOS dosimeter accessible again for manual changes.

WARNING

The XRS probe can be mechanically damaged or become unusable when the XRS X-ray source is inserted into the verification devices. This may also lead to incorrect irradiation.

- Always use the V-block guide when inserting the XRS into the verification devices.


WARNING

The XRS probe can be mechanically damaged or become unusable when the orientation of the V-block is incorrect during insertion of the XRS into the verification devices. This may also lead to incorrect irradiation.

- One side of the V-block is designed specifically to hold the XRS X-ray source, while the other side holds the verification devices. This is marked by an appropriate engraving on the V-block.


Key for test setup

- 1 IC shield
with IC ionization chamber (not visible in the illustration)
- 2 PAICH
with RTD resistance temperature device
- 3 XRS X-ray source
with IRM internal radiation monitor
- 4 X-block
- 5 UNIDOS E dosimeter
- 6 PRS500 control console

(Blank page, for your notes . . .)

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. A vertical red margin line runs down the right side of the page. The paper appears to be from a notebook or a standard sheet of stationery. There are no markings, text, or drawings on the page.

The test is performed as follows:

Activating the function

- In the "System Quality Assurance" mask, select the option "PAICH Output Check".

Video instruction

Press the  button to start an SQA or system calibration procedure, and the application software will visualize the steps to be performed in a video instruction with displayed text, see Page 136. The application software checks internally whether all verification devices required have been connected and are in position. If this is not the case, the procedure cannot be continued.

In the Dosimeter field, the color code for "Active" (3) switches to green if communication takes place between the UNIDOS E dosimeter and the user terminal

- If the UNIDOS E dosimeter has not yet been reset, you now have to click on <Zero> button (4) in the "Dosimeter" field.
- Click on <Start> button (11) to start the procedure.

There is an initial warm-up phase before the values are read out and averaged.

- Check that the IC is properly connected to the dosimeter, and that the UNIDOS E dosimeter is correctly set up, i.e. measuring unit set to pA, low range.

Refer to the UNIDOS E dosimeter user manual for more information.

The value indicated in "Dosimeter Reading [pA]" box (7) is the pure reading of the UNIDOS E dosimeter.


The value indicated in "Output Dose Rate [Gy/min]" (8) is the value computed internally with temperature and air pressure taken into account. The displayed numeric value is the UNIDOS E dosimeter reading averaged over time.

- Press <Accept> button (10) to confirm the value.

After accepting the measured value, you can start treatment planning (see Page 163) or perform a further verification step.

Overview of the "PAICH Output Check" screen mask

- 1 Selecting the "PAICH Output Check" verification step.
- 2 Display of temperature and air pressure
- 3 "Dosimeter, Active" color code display (communication activity)
- 4 <Zero> button, resetting the UNIDOS E dosimeter
- 5 Text display of the current step in the procedure
- 6 <Stop> button for stopping the verification procedure
- 7 "Dosimeter Reading" display, readings of the dosimeter UNIDOS E
- 8 "Output Dose Rate" display, computed value
- 9 <Reject> button for discarding the new values
- 10 <Accept> button for confirming the new values
- 11 <Start> button for starting the verification procedure
- 12 Graphic display of the beam deviation from the center

- Check the XY diagram on the bottom left of the screen mask which indicates the electron beam runout in the X and Y directions. The red circle represents the absolute limit value, i.e. the system is shut down when this value is exceeded.

In practical use, no more than 50% of this tolerance should be utilized. Therefore, you should recenter the XRS probe when approx. 50% of the admissible maximum runout are reached. Activate the "Probe Adjuster" function for this purpose, see Page 138.

The XRS probe can be mechanically damaged when the XRS is removed from the verification device.


- Always use the V-block guide when removing the XRS from the verification device.


"IRM Rate History" window, subsequent window to "PDA Source Check"

The diagram shows the time variance of the dose rate and the percentage deviation from the factory calibration.

- Press the <Yes> button to save the new values.
- 1 <Yes> button
The values are saved.
 - 2 <No> button
The values are discarded. You must repeat the system verification procedure (SQA).
 - 3 Display of percent deviation
Percent deviation of the current measurement from the first value ((Day 0)

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | Review Treatment | Calibration | 

Current report for the "PAICH Output Check"

After system verification using the "PAICH Output Check" routine, you can display or print out the current result, i.e. the latest SQA report of the "PAICH Output Check" routine:

Displaying reports. <Open Latest>

In the <Reports> menu, click on <Open Latest> and select <PAICH>. The report of the latest "PAICH Output Check" is displayed as a table in a separate window.

Exporting reports of the "PAICH Output Check", see Page 108.

The <Print> button permits you to print out the table as a PDF file and export it.

Use the <Done> button to close the window.

SQA PAICH Output Check	
Institution:	Carl Zeiss
Date/Time:	6/18/2007 14:43:58
Responsible Operator:	bess
XRS Serial Number:	000170
Beam Voltage [kV]:	50.0
Beam Current [uA]:	40.0
Dosimeter Model:	UNIDOS-E
Dosimeter Serial Number:	00250
Dosimeter Calibration Factor:	1.000000E+0
Ion Chamber Model:	PTW23342
Ion Chamber Serial Number:	N23342-1006
Ion Chamber Calibration Factor [Gy/C]:	1.141000E+9
PAICH Model:	PAICH450
PAICH Serial Number:	507190
Pressure [kPa]:	95.5
Standard Pressure [kPa]:	101.3
Temperature [°C]:	21.4
Standard Temperature [°C]:	22.0
P&T Correction Factor:	1.058
Output Value [Gy/min]:	2.873
Dose Rate Correction[%]:	6.0
Result:	Pass

Only for predecessor systems: checking the ERM external radiation monitor

Only predecessor systems enable you to check the function of the ERM external radiation monitor. The count rate of the ERM external radiation monitor is checked against the count rate of the IRM internal radiation monitor. The ERM count rate obtained during this verification procedure will not be used in treatment administration.

The "ERM Check" procedure is not enabled until the "PDA Source Check" has been successfully completed, i.e. the XRS probe must be straightened in accordance with the relevant specifications prior to the ERM check. The isotropy of the X-ray field must have been adjusted to be within the required limits.

This image shows a single page of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There is no handwriting or other markings on the paper.

Entering treatment parameters

"Enter Prescription" screen mask

The "Enter Prescription" screen mask permits you to select the components of your equipment to be used for the current treatment session, e.g. an XRS X-ray source.

Activating the function

- Click on the "Enter Prescription" tab.

Filling in the screen mask

- Select the XRS X-ray source in the "XRS No." box.
The XRS must have been successfully verified to be selectable here, see "Equipment Inventory" screen mask.
- Check the XRS No. to make sure that you have selected the right XRS.
- Select the voltage to be used for treatment in the "XRS Voltage" box.
You can only select a voltage for which the selected XRS X-ray source has been successfully verified, see "Equipment Inventory" screen mask.
- Select the current to be used for treatment in the "XRS Current" box.
You can only select a current for which the selected XRS X-ray source has been successfully verified, see "Calibration"/"IRM Linearity" screen mask.

Select the treatment method by setting either "Spherical" or "None" for the "Applicator Type".


For tumor bed irradiation, you require an applicator; therefore select "Spherical". For interstitial treatment, no applicator is required; therefore select "None". (Only applicators for 50 kV are currently available.)

- If you need an applicator for a treatment session and have selected "Spherical" for the Applicator Type, you have to set the applicator size in the "Size [cm]" box:
You can only select applicators which have been set to "Active" in the "Equipment Inventory" screen mask. If "Equipment Inventory" contains several applicators with the same size, you have to select the applicator to be used on the basis of its serial number (Applicator No. box).
- Check the serial number to make sure that you have selected the right applicator.
- Enter the patient information by clicking on the "Patient ID" box. Use the <More> button to open a follow-up window for entering patient data, see Page 168.

**WARNING**

If the user confuses the treatment plans, radiotherapy will be performed with incorrect treatment data. This can lead to injury through radiation.

- Make sure that you have selected the right treatment plan.

- Enter the dose to be administered in the target area in "Prescription Dose" box (18). The dose in the target area results from the distance from the isocenter or the distance of the applicator surface from the target area ("Treatment Depth" box (17)).

About the screen mask

The lower part of this screen mask shows the dose/depth curve which indicates the dose present at specific depths (mm). It permits you to check to which extent adjacent tissue is irradiated.

Two sets of crosshairs are displayed in the diagram.

- Crosshairs (16) are stationary in the position of the planned treatment parameter. For the values, see boxes (17) and (18).
- The second crosshairs (15) can be moved using the mouse. The values of the second crosshairs (15) are displayed at the top right in two boxes (13).

In "Depth [mm]" box (13), you can also enter a precise depth value. The second set of crosshairs (15) will then jump to this value in the diagram.

After you have entered the dose to be administered at the relevant depth, the treatment time is displayed. Check the treatment time in view of the patient's tolerance of the anesthetic.

After successful entry of all treatment parameters required, color code "Green" is displayed for "Planning Complete" (12). You can now check the treatment parameters using the "Verify Treatment Parameters" procedure (Page 173).

"Enter Prescription" screen mask

- 1 Selecting the XRS X-ray source
The selection is made manually via the serial number.
- 2 Selecting the accelerating voltage
- 3 Selecting the beam current
Only currents producing a minimal count rate in the determination of IRM linearity (see) are entered in the calibration file and can be used for treatment planning.
- 4 Selecting the type of treatment
- 5 Selecting the applicator size
- 6 Serial number of the applicator
If your applicator set (Inventory) contains several applicators with the selected size, you have to select the serial number of the applicator to be used.
If your applicator set contains only one applicator with the selected size, its serial number is entered automatically.
- 7 Continue with previous patient
Tick the checkbox if treatment of the same patient is continued without system reverification within 90 minutes of the first treatment.
- 8 Patient ID, entry and display of the patient ID
The file name of the treatment record is automatically generated from the patient ID entered here and the current date. If you wish to use an additional search option, e.g. via the patient's name, you have to enter the search item (e.g. the patient's name) in this box, after the Patient ID.
- 9 <More> button
Use the <More> button to open the follow-up window "Patient Information" for the entry of further patient data, see Page 168.
- 10 Treatment number, "Treatment #" display box
The INTRABEAM System automatically allocates a number for each treatment. The numbers are unique and cannot be changed.
- 11 Internally computed treatment duration and dose rate
- 12 "Planning Complete" color code display,
lights up green after all treatment parameters required have been entered.
- 13 Planned treatment parameters (depth/dose) are displayed
You can enter the required depth in the depth field.

Equipment Inventory | System Quality Assurance | **Enter Prescription** | Verify Treatment Parameters | Treatment | Review Treatment | Calibration

14 Graphic display of the depth dose rate

15 Crosshairs, movable with the mouse function


In the depth dose rate diagram, crosshairs are displayed at the planned dose / depth treatment parameter. You can use the mouse to shift these crosshairs along the depth / dose curve.

The values at the crosshair position (depth/dose) are displayed at the top right of the diagram (13). You can enter the required depth in the depth field. The values displayed indicate which dose is administered at which depth in the tissue (17) and (18) during the planned treatment.

16 Fixed point in accordance with the treatment parameters (17 and 18)

17 Entry of the treatment depth, "Treatment Depth" box

18 Entry of the dose to be administered, "Prescription Dose" box

Follow-up window "Patient Information"

The "Patient Information" window is used to enter patient data for the envisaged treatment.

Activating the "Patient Information" window

- Press the <More> button in the "Enter Prescription" screen mask, see Page 167, Pos. (9).

Entering patient data for the current treatment session

- Fill in all input boxes of the Patient Information window.
- Note the input specifications for text entries:
 - Do not use any special characters such as < > | * " ? : / \
 - Do not enter more than 20 characters in a text box.

"Patient ID" and "Treatment #" boxes

The information in the "Patient ID" and "Treatment #" boxes cannot be edited. The Patient ID and Treatment No. were defined in the "Enter Prescription" screen mask, Page 167.

"Application" menu for detailed application information

"Application" (1) offers you a pull-down menu for the selection of different applications:

None
Left Brain
Right Brain
Left Breast
Right Breast

An overview image for the respective application is shown in display window (2).

Overview of "Patient Information" follow-up window

- 1 "Application" menu
Selection of detailed information on the different applications.
- 2 Display field
Display of an overview image for the application selected in "Application" menu (1).
- 3 Display field
Detailed image for marking the tumor position
- 4 <OK> button
Press the <OK> button to confirm that you wish to use the patient data currently entered.
- 5 <Cancel> button
Press the <Cancel> button to discard the patient data currently entered.

The screenshot shows the 'Patient Information' window with the following fields and controls:

- 1** Application: A dropdown menu currently set to 'None'.
- 2** Display field: A large empty rectangular box for an overview image.
- 3** Display field: A detailed coordinate grid for marking tumor position. The grid has X and Y axes labeled 'Conditional Units' ranging from 0 to 10. Above the grid, 'Cursor Position' is shown with X=0.0 and Y=0.0, and a 'Lock/Unlock' checkbox.
- 4** Done button: A button at the bottom left.
- 5** Cancel button: A button at the bottom right.

Other visible fields include:

- Patient ID: ID1234
- Treatment #: 404
- Other ID
- Given Name
- Family Name
- Insurance ID
- Sex: Female
- Date of Birth (MM/DD/YYYY)

At the bottom, there is a 'Use New Values?' checkbox.

Entering detailed patient information

Display window (3) of the Patient Information screen mask permits you to enter the tumor position (coordinates) and to branch to other detailed patient information.

Activation

Use the "Application" menu (Page 168, item 1) to select the detailed patient information required. The information is shown in the display field (Page 168, item 2).

The example "Left Breast" with detailed patient information is explained in the following.

Entering the tumor position (coordinates)

Using the cursor function, you can mark the tumor position in the drawing. Cursor (4) can be moved using the mouse if "Lock/Unlock" (2) has been set to "Unlock", i.e. the checkbox has not been ticked. Untick the checkbox, if necessary.

- Place cursor (4) at the tumor position. The X and Y coordinates are displayed above the drawing (3).
- Lock the cursor by setting "Lock/Unlock" (2) to "Lock". Tick the checkbox for this purpose.

The "Lock" setting prevents inadvertent shifting of the cursor.

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
Even if "Continue with previous patient" (Pos. 6) has been selected in the "Enter Prescription" screen mask (ticked checkbox), you can edit the entries in the following boxes: "Other ID", "Treatment #" and "Cursor Position".

Overview of the Patient Information window

- 1 "Application" menu
You can switch to different detailed patient information
- 2 "Lock/Unlock"
You can lock the cursor in position (tick the checkbox) to avoid inadvertent shifting of the cursor.
- 3 X/Y display
The XY coordinates of cursor (4) are displayed.
- 4 Cursor function
Used to enter the tumor position (coordinates). The cursor can be moved using the mouse.
- 5 <Cancel> button
Discarding the current entries.
- 6 <OK> button
Using and saving the current entries.

This image shows a single page of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page. There is no handwriting or printed text on the page.

Checking the treatment parameters

"Verify Treatment Parameters" screen mask

The "Verify Treatment Parameters" screen mask is used by the radiotherapist (treating physician) to check and confirm the treatment parameters entered by the physicist in the preceding "Enter Prescription" screen mask.

The radiotherapist (treating physician) acknowledges the correctness of the entered treatment parameters by the entry of his/her personal password.

Data can no longer be edited in the "Verify Treatment Parameters" screen mask. If any data is to be changed, you have to return to the "Enter Prescription" screen mask.

The "Verify Treatment Parameters" screen mask can only be activated if the preceding "Enter Prescription" screen mask has been completely filled in and the color code display for "Prescription Complete" is green.

Activating the function

- Activate the "Verify Treatment Parameters" screen mask by clicking on the <Verify Treatment Parameters> tab.

Password Entry

The entry of the password is equivalent to a signature. Depending on local legislation, both the radiotherapist (treating physician) and the physicist may be required to personally sign the treatment record which is output by the INTRABEAM system at the end of the treatment session.

Password regulation

- If the physicist logs into the system, this is equivalent to a signature.
- Even if the physician has logged into the system first, the following applies:
the physicist has to sign, followed by the physician.

Only if the verification of the treatment parameters has been correctly completed will you be able to activate the "Treatment" screen mask and conduct the treatment.

Overview of the "Verify Treatment Parameters" screen mask

- 1 **Checking the displayed treatment parameters**
The treatment parameters entered by the physicist in the "Enter Prescription" screen mask are displayed here and must be checked and acknowledged by the radio oncologist (treating physician).
- 2 **User log-in: physician and physicist**
The radio oncologist (treating physician) acknowledges the correctness of the entered treatment parameters by the entry of his/her personal password.
If several persons have been entered in the system, you can select your name from a pull-down menu.
- 3 **<Verify> button**
By pressing the <Verify> button, the radio oncologist (treating physician) confirms that the entered treatment parameters are correct and closes the verification of the treatment parameters.
- 4 **"Verification Complete" color code display**
After the radio oncologist (treating physician) has pressed the <Verify> button, the color code display lights up green, indicating that the verification of the treatment parameters has been correctly completed. Only then will it be possible to activate the subsequent "Treatment" screen mask.
- 5 **<More> button**
Use the <More> button to open a follow-up window for the display of patient data, see Page 168.

Equipment Inventory | System Quality Assurance | Enter Prescription | **Verify Treatment Parameters** | Treatment | Review Treatment | Calibration

1 {

XRS		Applicator		Patient Information	
No. 000170	Voltage 50 [kV]	Type Spherical		Patient ID	ID1234
	Current 40 [μA]	Size[cm] 3.0			More
XRS Name	1st XRS 50	No. 265401		Treatment #	404

2 {

Prescription Dose	20.00 [Gy]	Treatment Time	0:23:49 [h:m:s]
Treatment Depth	0.00 mm from Surface	Dose Rate at Depth	0.840 [Gy/min]

Prescribing Physician	Password	Ok	Responsible Physicist
Physician	*****		Physicist

3 {

4 {

5 {

This image shows a single page of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page. There is no handwriting or printed text on the page.

Operation

Checklist for INTRABEAM Core System	178
Conducting treatment	180
Troubleshooting	191
Documenting the number of applicator sterilizations	196
Review treatment, displaying treatment parameters	199

Checklist for INTRABEAM Core System

- Always check the following points before using the instrument(without patient!):

Readiness for operation

- ✓ The components were properly washed, disinfected and sterilized, if required.
- ✓ The INTRABEAM Core System with the XRS X-ray source has been verified.
- ✓ Calibration is still valid.

Accessories

- ✓ The other equipment, e.g. applicators, sterile drapes (INTRABEAM Drapes) etc., have been examined and found to be complete sets using the corresponding Instructions for Use.
- ✓ The use-by-date has not expired.
- ✓ The calibration data of the applicators has been entered into the application software.

INTRABEAM Cart (equipment cart)

The INTRABEAM cart is

- ✓ equipped with the components and accessories required for treatment,
- ✓ properly positioned, i.e. it is standing at least 1.5 m away from the OR table (patient),
- ✓ secured in position with locking tabs on the front casters, preventing it from rolling away.

You will find the "Equipment checklist for the operating room" in the chapter "Supplementary documents", on Page 265.

(Blank page, for your notes . . .)

Conducting treatment

**WARNING**

If X-ray generation starts before the XRS probe is inserted in the target area, unshielded radiation is emitted that injures the patient and user.

- Do not activate the XRS before the X-ray probe with applicator is in the target area.

**WARNING**

If critical organs are not sufficiently shielded during the treatment, they may be injured through X-radiation.

- Please use shields, to screen off critical organs.

**WARNING**

If the patient is not sufficiently shielded against radiation or if no adequate measures have been taken to protect the patient from radiation, X-radiation may be emitted into the surrounding area and injure the user.

- Please ensure the target area is shielded off suitably as a result of using shields.
- If required, use portable radiation shields to protect OR staff against radiation.
- Stay away from the target area as far as possible.

**CAUTION**

The probe tip of the XRS X-ray source contains non-compliant material. If the user applies the x-ray probe directly into the tissue, incompatible material will have contact with the wound tissue, which can result in allergic reactions.

- Never insert the XRS probe directly in the tissue.
- Always mount an applicator on the XRS probe.

Radiating a tumor bed

- Mount the XRS X-ray source with the attached probe guard on the stand.
- Connect all the required cables on the INTRABEAM system, see Page 92. Connect the cables of the XRS x-ray source on the INTRABEAM stand with the control unit.
- Mount the applicator selected in treatment planning on the stand.

CAUTION

Bacterial contamination in the sterile OR area can occur if a drape is not used.

- In general, always use a drape on the INTRABEAM stand in the OR.

- Cover the INTRABEAM stand with a sterile drape, see Instructions for Use "NC32 INTRABEAM floor stand"..

"Treatment" screen mask

Use the "Treatment" screen mask to conduct treatment.

Before activating the "Treatment" screen mask, make sure that the system has been set up completely and, in particular, that all cable connections have been made.

The "Treatment" screen mask can only be activated if the "Verify Treatment Parameters" mask has been completely filled in and the color code display for "Planning Complete" is green, see Page 173.

Activating the function


- Activate the "Treatment" screen mask by clicking the "Treatment" tab.
- Make sure the XRS X-ray source has been positioned as required for the type of treatment involved.

Overview of the "Treatment" screen mask

- 1 Interlock Status
The interlock status indicates whether all interlocks of the system have been closed, or the current status of the interlocks.
- 2 IRM/ERM diagram
Current values (expected/actual) of the IRM and ERM (of predecessor systems only), if connected, are displayed.
- 3 Interruption in the diagram
The diagram shows an interruption in treatment by <Pause> button (4) and continuation of treatment by <Resume> button (4).
- 4 Ready / Start / Pause / Resume button
The labeling of button (4) changes depending on the current mode of operation:

Ready	In the standby mode during system start.
Start	Start phase completed; system is ready to start radiation emission.
Pause	Interrupts radiation emission while treatment is in progress.
Resume	Continuation of treatment after a pause.
- 5 Stop button for aborting treatment
- 6 Opening the input window
Window for the entry of comments during treatment; for details see Page 186.
- 7 Display window
Text display of interactions and system messages
- 8 Graphic display of beam deflection
- 9 Percentage deviation from the expected IRM rate
- 10 For predecessor systems only:
Values of ERM external radiation monitor
Value = 0 means: ERM not connected

- 11** Current IRM value
- 12** For predecessor systems only:
Values of ERM external radiation monitor
Value = 0 means: ERM not connected
- 13** Treatment #
Number of the current treatment session, automatically assigned by the INTRABEAM system

Starting treatment

The INTRABEAM system checks all interlocks and that all cables and the X-ray source have been connected.

If these requirements are not met, button (4) is disabled.

- If all requirements for treatment have been met, press the <Ready> button (4).

The user terminal transmits the treatment data to the PRS 500 control console. As soon as the PRS 500 control console is ready, the labeling of button (4) changes from <Ready> to <Start>.

You must no longer make any changes to the system at this stage, such as unplugging any cable connections.

- To activate radiation emission: press <Start> button (4).
Once treatment has been started, the labeling of button (4) changes from <Start> to <Pause>.
- **The following data must be monitored during treatment:**

The current IRM values should correspond to the set values within a tolerance of $\pm 10\%$. If the current values are out of tolerance, the application software displays an error message and stops treatment.

- Observe the beam deflection in "Beam position" window (8). If the current values are out of tolerance, irradiation is stopped and an error message is displayed.

In most cases, the current beam position moves slowly to the edge of the circle. This is frequently due to tensioning of the XRS probe caused by patient or tissue movement. A brief pause to relieve the tension in the XRS probe may help to prevent stoppage of the treatment.

If an ERM external radiation monitor is used, the activation of the <Start> button and every activation of the <Resume> button initiates the determination of the ERM count rate averaged over time in proportion to the count rate of the IRM internal radiation monitor, with the current count rate of the ERM external radiation monitor being 100%.

Overview of the "Treatment" screen mask

2 Diagram

Current values (expected/actual) of the IRM and ERM, if connected, are displayed.

3 Drop in the IRM rate, caused by a pause

4 Ready / Start / Pause / Resume button

The labeling of button (4) changes depending on the current mode of operation:


Ready In the standby mode during system start.

Start Start phase completed; system can start radiation emission.

Pause Interrupts radiation emission

Resume Continuation of treatment after a pause.

8 Graphic display of beam deflection

Adding comments during treatment

You can enter comments of any type in the "Treatment" screen mask during treatment.

Opening the input window

- Click on quill icon (12).
Input window (11) is opened.

The open input window (11) is superimposed on the display window (display window see Page 183, item 7).

You can enter your comments on the treatment in progress in the open input window (11).

You can continue to enter comments in the separate input window even if the treatment session has been interrupted (<Pause>) or stopped (<Stop>).

Closing the input window

Input window (11) is minimized as soon as you change to a different field in the "Treatment" screen mask, or:

- When input window (11) is open, you can click on the keyboard symbol to minimize the window.

The display window (see Page 183, item 7) is visible again.

Saving comments

Your comments are saved automatically. They are displayed in the "Review Treatment" screen mask under "Comments", but they can no longer be edited there.


Overview of the "Treatment" screen mask with the input window

11 Input window

You can use this window to enter comments of any type during treatment.

12 Quill icon

Click on this icon to minimize input window (11).

11

12

Interrupting treatment (Pause mode)

A pause in the treatment session may be required to make minor adjustments, e.g. to correct the patient's position, or if a beam deflection was detected in the "Beam Position" window (8). A pause in the treatment session also permits you to slightly correct the position of the XRS X-ray source in order to reduce any force effect on the XRS probe. A brief pause to relieve the tension in the XRS probe may help to prevent stoppage of the treatment.

- Use <Pause> button (4) to temporarily interrupt X-radiation emission. A long beep is heard.
Button (4) is labeled <Resume> during the pause.
- Press <Resume> button (4) to continue treatment. When treatment is resumed, button (4) is labeled <Pause> again.
The ERM external radiation monitor (only concerns predecessor systems) is re-initialized after a pause, with the previous values being appropriately scaled.

The control panel of the PRS 500 control console also includes a <PAUSE> button which permits you to interrupt the treatment session. However, you cannot continue treatment by once again pressing the <PAUSE> button on the PRS 500 control console; instead, you have to press the <Resume> button (4) on the user terminal.

Stopping treatment

Press <Stop> button (5) to abort the treatment session. The INTRABEAM system indicates the dose administered up to this point and the difference from the nominal dose in percent.

If you abort treatment using <Stop> button (5) or the Emergency Stop switch, you must re-verify the XRS X-ray source, unless you continue to treat the same patient using the "Continue with previous patient" option. In that case, you can administer the remaining dose within 90 minutes without having to re-verify the XRS X-ray source, see the "Enter Prescription" screen mask on Page 167

Overview of the "Treatment" screen mask

4 Ready / Start / Pause / Resume button


The labeling of button (4) changes depending on the current mode of operation:

Ready	In the standby mode during system start.
Start	Start phase completed; system can start radiation emission.
Pause	Interrupts radiation emission while treatment is in progress.
Resume	Continuation of treatment after a pause.

5 Stop button

for aborting treatment

8 Graphic display of beam deflection

End of the treatment session

Treatment has been completed as planned when the administered radiation dose is identical to the prescribed radiation dose (see Treatment planning). If the administered dose deviates from the prescribed value (e.g. if treatment was stopped), you have to plan a further treatment session with new parameters to administer the remaining dose.

Automatic backup after treatment

All patient files, verification reports and treatment reports are stored in a database under D:\Databases.

The following files can be accessed under D:\Databases:

- ZEISS PRS500 EI.mdb
 - comprising data on equipment and verifications of the XRS
- ZEISS PRS500 Patients.mdb
 - comprising patient and treatment data

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If you wish to treat the same patient again, you can plan a further treatment session within 90 minutes without having to reverify the system.

- Tick the checkbox for "Continue with previous patient" when entering the treatment parameters in the "Enter Prescription" screen mask, see Page 166.

Troubleshooting

NOTE

After a breakdown, the User may not be able to restart the INTRABEAM core system without information from the Instructions for Use and will have to discontinue the treatment.

- Please ensure that the user is specially trained on the subject of "Troubleshooting" and that a copy of the Instructions for Use are always handy.

Failure of the user terminal

If the connection between the control unit and the user terminal is interrupted after a treatment has already been started, e.g. as a result of detaching the connection cable, interrupting the power supply to the user terminal, or similar, then the control unit will continue the treatment on its own.


If the user discontinues the treatment due to a breakdown in the user terminal, then the treatment time will extend itself unnecessarily.

- Allow the INTRABEAM system to continue, do not under any circumstances restart the application software.

If the control unit continues the treatment on its own, then all the information required for the treatment will appear on the operating panel display:

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The PRS 500 control console will end the treatment session after the prescribed dose has been administered.

- Restart the application software of the user terminal after the PRS 500 control console has ended the treatment.
The PRS 500 control console and the user terminal automatically exchange the treatment data so that all data can be checked and output on the user terminal in the "Review Treatment" module of the application software.

Failure of the PRS500 control console

If the power supply to the PRS 500 control console is interrupted during a treatment session or if the device is switched off by pressure of the Emergency Stop switch, the treatment data existing up to this point is automatically saved. Proceed as follows in this case:

- End the application software by selecting <File> / Exit (proper shutdown of the application software)
- Switch the PRS 500 control console back on.
- Wait until the self-test is completed and restart the application software.

The PRS 500 control console and the user terminal automatically exchange the treatment data so that all important data can be checked and output on the user terminal in the "Review Treatment" module of the application software.

In this case, less than the planned dose has been administered. The remaining dose can be administered in a further treatment session with new parameters. If the remaining dose is to be administered with the INTRABEAM system, the XRS must be re-verified. Alternatively, the remaining dose can be administered without reverification by using the "Continue with previous patient" function.

Automatic pause mode

If an interlock error occurs during treatment, e.g. because

- the applicator is no longer in the INTRABEAM interface,
- the external interlock (door contact) is interrupted,
- the XRS cable has been disconnected
- or due to another reason

the INTRABEAM system automatically goes into the pause mode.

You can see the current interlock status at the top left of the "Treatment" screen mask.

- Properly connect all interlocks until the interlock status display "TREATMENT" appears.
- Then press the <Resume> button and continue treatment.

Troubleshooting in the event of ABC errors

(ABC-Error = Automatic Beam Centering Error)

An error message may be displayed if the XRS probe is excessively deflected. The cause is frequently a twisting of the x-ray probe, caused by a position change in the patient or the tissue.

If the current values are outside of tolerance, the radiation treatment will be discontinued and an error message will be shown on the monitor.

If there is mechanical tension on the applicator, you can watch this during the radiation on the graphic depiction of the beam deflection (on the screen mask "Treatment", screen "Beam position") and ascertain that the current beam position is gradually wandering towards the edge of the circle. In this case, stop the treatment session by pressing the <Pause> button. Rectify the mechanical tension on the applicator by briefly detaching the brake on the INTRABEAM floor stand.

This requires a sterile staff member!

- When releasing the brake, pay attention to the possible movement of the stand arm.

If an error occurs:

- Note the prescribed dose and the dose administered up to this point.
- Determine the dose still to be administered.
- Release the stand brakes using the hand switch to relieve the tension in the XRS probe.
- Switch to the "Enter data" tab and use "Continue same patient" to plan this differential dose as a new treatment and continue with the radiation.

If the error occurs again:

A sterile staff member is required (sterile person) to remove the XRS x-ray source from the treatment field.

- Switch tables in "System Quality Assurance (SQA)."
- Align the XRS tip and carry out the other routines to check the system (SQA): Dynamic Offset, PDA Source Check and PAICH Output Check.
- Plan a new treatment session for the remaining dose.

If the error persists:

- Notify Carl Zeiss Service.

Treatment stopped for other reasons

- Note the prescribed dose and the dose administered up to this point.
- Determine the dose still to be administered.
- If required, restart the INTRABEAM System again.
- Switch to the "Enter data" tab and use "Continue same patient" to plan this differential dose as a new treatment and continue with the radiation.

(Blank page, for your notes . . .)

Documenting the number of applicator sterilizations

Only the applicators selected in "Enter Prescription" are covered by the system's automatic incrementation function once the treatment has been started. The sterilization (incl. cleaning and disinfection) of the other applicators has to be documented manually using the menu item <Applicator Sterilization Monitor>.

Several applicators are often provided in the OR for radiation therapy treatment. All these applicators are then considered to be "contaminated" and need to be resterilized.


The "Applicator Sterilization Monitor" is a tool for monitoring the sterilized applicators. However, this tool does not release you from the responsibility of documenting the number of sterilizations for all applicators.

WARNING

Exceeding the maximum admissible number of sterilizations can damage the applicator and lead to incorrect radiation.

- Make sure that all sterilization cycles are counted and that no spherical applicator is used that has been sterilized more than 100 times.

Activating the function

Click on "Tools" in the menu bar and select the menu item "Applicator Sterilization Monitor". The "Applicators" window is displayed.

Selecting applicators

The "Applicators" window contains two lists. The "Available" list on the left includes the active applicators in the system, the "To Sterilizer" list on the right the applicators to be sterilized.

- Click on the applicator in the "Available" list that has to be sterilized.
- Move the selected applicator from the left to the right list by pressing the arrow button ">".

You can select several applicators at a time if you hold down the <Ctrl> button during selection. All selected applicators can then be moved to the right list using the arrow button ">".

Arrow button "<" permits you to correct your selection.

- When you have selected all "contaminated" applicators and moved them to the right list, click on the <Done> button once.

When the <Done> button is pressed, the system increments the number of sterilizations for the relevant applicators.

All applicators of your INTRABEAM system are managed in the "Applicator In-


ventory" table of the "Equipment Inventory" screen mask, see Page 112. As soon as one of the spherical applicators has reached the limit of 100 sterilizations, the relevant applicator is set to "inactive" in the "Equipment Inventory" screen mask and can no longer be used.

(If you wish to continue using a spherical applicator beyond the recommended maximum of 100 sterilizations, you can reactivate it. Please note, however, that you have to confirm this with your password, see Page 128.)

The following definitions have been made to avoid any double counts:

- Selected applicators are marked with a tick and cannot be selected again after reactivation of the "Applicators" window; they are blocked.
- Blocked applicators can be selected again after 60 minutes and after new login.

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Review treatment, displaying treatment parameters

"Review Treatment" screen mask

The "Review Treatment" screen mask permits you to review the course of treatment and to print out the treatment record.

The data of the treatment session currently completed is automatically saved in a file whose name is automatically generated from the data entered in the "Patient ID" box. For further details, see Page 166.

The "Review Treatment" screen mask automatically gives an overview of the data of the treatment currently completed.

Activating the function

- Activate the "Review Treatment" screen mask by clicking on the <Review Treatment> tab.

You can also display saved records from earlier treatment sessions. Use the <Load File> button to load a saved treatment record to the "Review Treatment" screen mask.

The <View File> button permits you to view the treatment record currently loaded.

Use the <Print Short Form Report> button to print out a brief record. The chapter "Data export" on Page 108 explains how to export treatment records from your user terminal.

Overview of the Review Treatment screen

- 1 <Report> button
Displays a detailed view (report) of the treatment data currently loaded in a separate window. For an example, see Page 204.
- 2 <Load data> button
To load a saved treatment report of an earlier treatment into the screen mask "Review".
The <View File> button then permits you to view the complete treatment record.
- 3 <More> button
The <More> button is used to open a subsequent screen to display patient data, see Page 168.
- 4 Display of the planned treatment time "Treatment time [Plan]"
Value calculated prior to beginning of treatment based on a 100% IRM Rate.
- 5 Display of the actual treatment time "Treatment time [Act.]"
The actual treatment time without possible break periods is displayed for the treatment that was just loaded. For explanations on the "Pause" function, see Page 188.
If the dose rate is changed during treatment, the actual treatment time may differ from the estimated treatment time by a maximum of 10%.
- 6 Tab row re the display field (6)
Various information on the treatment that was just loaded can be displayed in the display field (6) by clicking on a tab. For additional explanations on the Tabs, see Page 203.
- 7 Display field
Depending on the selected tab (5), you can also have various information on the treatment that was just loaded displayed here.
The display window in the illustration shows the Treatment Event Record.
- 8 <Print Short Form Report> button
To print out the short version of the treatment report for the treatment that was just loaded. This short form report is two pages in length.

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | **Review Treatment** | Calibration

Zeiss INTRABEAM System
File Reports Tools Help

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | **Review Treatment** | Calibration

XRS		Applicator		Patient Information	
No. 000170	Voltage 50 [kV]	Type	Spherical	Patient ID	ID1234
	Current 40 [μA]	Size[cm]	3.0		More 
XRS Name	1st XRS 50	No.	265401	Treatment #	404
Load File 	Prescription Dose	20.00	[Gy]	Treatment Started	11:01:39
View File 	Dose Delivered	20.00	[Gy]	Estimated Treatment Time	0:23:49 [h:m:s]
	Treatment Depth	0.00	[mm]	Actual Treatment Time	0:23:50 [h:m:s]
1	Physicist	Physicist	Physician	Physician	
Treatment Event Record Depth-Dose IRM/ERM Beam Position Comments					
11/08/2007 11:01:19 Treatment screen selected Responsible Physicist: Physicist Prescribing Physician: Physician Patient ID: ID1234 Treatment #: 404 11:01:20 Checking interlock conditions, please wait ... 11:01:20 Filament check passed 11:01:21 All interlock conditions are satisfied. 11:01:26 Loading Treatment Parameters into the CC, please wait ... 11:01:29 Treatment Parameters loaded. 11:01:30 READY state reported by the CC 11:01:38 Operator pressed 'Start' button 11:01:39 Treatment started 11:06:38 Operator pressed 'Pause' button 11:06:38 Treatment paused 11:07:00 Operator pressed 'Stop' button					
				 Print Short Form Report 	

Overview of the tab bar for the display window

- 9** Tab bar above the display window
By clicking on the tabs in the tab bar above the display window, you can view different data from the treatment session currently loaded in the display window.
- 10** Display window
Different data from the treatment session currently loaded can be shown in the display window. Select the information required by clicking on the relevant tab in the tab bar above the display window.
- 11** <Treatment Event Record> tab
The course of treatment is displayed in the form of a list in the display window, as shown in the illustration.
- 12** <Depth Dose> tab
The depth-dose curves (depth [mm] / dose [Gy], see Page 167) of the planned and the actually administered dose are shown in the display window.
- 13** <IRM/ERM> tab
The display window shows the current values (expected/actual) of the IRM and ERM, if connected, see Page 182.
- 14** <Beam Position> tab
The display window shows the beam deflection over time.
- 15** <Comments> tab
The display window shows the comments manually entered during the treatment session.

Zeiss INTRAM System									
File Reports Tools Help									
Equipment Inventory		System Quality Assurance		Enter Prescription		Verify Treatment Parameters		Treatment	
						Review Treatment		Calibration	
XRS		Applicator				Patient Information			
No. 000170	Voltage 50 [kV]	Type	Spherical		Patient ID	ID1234			
	Current 40 [μA]	Size [cm]	3.0		More 				
XRS Name	1st XRS 50	No.	265401		Treatment #	404			
Load File 	Prescription Dose	20.00	[Gy]	Treatment Started	11:01:39				
	Dose Delivered	20.00	[Gy]	Estimated Treatment Time	0:23:49 [h:m:s]				
View File 	Treatment Depth	0.00	[mm]	Actual Treatment Time	0:23:50 [h:m:s]				
Physicist		Physicist		Physician			Physician		
9	Treatment Event Record	Depth-Dose	IRM/ERN	Beam Position	Comments				
10	11:08:00 11:01:19 Treatment screen selected Responsible Physicist: Physicist Prescribing Physician: Physicist Patient ID: ID1234 Treatment #: 404 11:01:20 Checking Interlock conditions, please wait 11:01:20 Filament check passed 11:01:21 All Interlock conditions are satisfied. 11:01:26 Loading Treatment parameters into the CC please wait ... 11:01:29 Treatment Parameters loaded. 11:01:30 READY state reported by the CC 11:01:38 Operator pressed Start button 11:01:39 Treatment started 11:06:38 Operator pressed Pause button 11:06:38 Treatment paused 11:07:40 Operator resumed treatment								
 Print Short Form Report 									
	11	12	13	14	15				

Viewing and printing the treatment record

When you press the <View File> button in the "Review Treatment" screen mask, the detailed treatment record of the current treatment session or of the treatment session currently loaded is displayed in a separate window.

Creating and saving PDF files

Select the "CutePDF Writer" printer in the printer selection. This "printer" creates a PDF file. The "SAVE" function permits you to save the PDF file to your system.

Recommended storage location: D:\TreatmentRecords

Use the <Print> button to print out the record, see the printout example on the opposite page.

Use the <Done> button to close the window.

Example of a treatment record, page 1

[illegible]


Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | Review Treatment | Calibration

Treatment record, printout example


INTRABEAM System ZEISS

Treatment Record

Institution	Carl Zeiss
Patient Given Name	Given Name
Patient Family Name	Family Name
Patient Date of Birth	16-Nov-1974
Patient Sex	Female
Patient ID	ID1234
Treatment #	404
Insurance ID	4567A
Other ID	7654B
Treatment Date	08-Nov-07
Treatment Start/Stop Time	11:01:39 / 11:26:33
Prescribing Physician	Physician
Responsible Physicist	Physicist
XRS Serial Number	000170
XRS Name	1st XRS 50
Voltage [kV]	50.0
Current [μA]	40.0
Applicator Type	Spherical
Applicator Size [cm]	3.0
Applicator Serial Number	265401
Applicator Sterilizations to Date	4
Applicator Sterilizations Max	100
Estimated Treatment Time	0:23:49
Actual Treatment Time	0:23:50
Prescription Dose [Gy]	20.00
Prescription Depth [mm]	0.00
Dose Delivered [Gy]	20.00
Deviation from Prescription [%]	0.0
PDA Report Name	PDA_XRS000170@50.0kV_200711081052
PAICH Report Name	PAICH_XRS000170@50.0kV_200711081052

Application Type: Left Breast

[illegible]

Care and maintenance

Calibration regulations	209
Calibration	210
Date and time synchronization	212
Determining the deflection sensitivity	214
Determining the linearity of the IRM internal radiation monitor	218
Calibrating the air pressure sensor	220
Calibrating the temperature sensor	223
For predecessor systems only: Dilator Adjustment.....	226
Troubleshooting table	227
Changing the line fuse.....	233
Resterilization.....	234
Cleaning and disinfecting components	235
Cleaning the touchscreen	237
Cleaning the equipment trolley(INTRABEAM Cart)	238
Disinfecting user terminal and cart.....	238
Handling, transport and shipment	239
Disposal.....	240
Ordering data	241
Spare parts	242
Accessories	243

[illegible]

Calibration regulations

The following calibration procedures are vital for safe, reliable operation of the INTRABEAM System:

Calibration of XRS X-ray sources

In the calibration process, the depth-dose rate characteristics of the XRS radiation field are quantified, i.e. the XRS output and radiation characteristics are calibrated.

WARNING:

Non-compliance with the calibration intervals for the XRS X-ray radiation source can lead to incorrect radiation.

- Have the XRS X-ray radiation source calibrated every year (recommended by the manufacturer ZEISS). If shorter calibration intervals are stipulated by national codes and regulations, you must comply with them.
- Have your XRS X-ray radiation sources calibrated in the ZEISS manufacturing plant. As part of this calibration, a new calibration file is created for the XRS and delivered to the user along with the XRS.

The calibration interval starts with the initial commissioning. As there can be a period of several months between the calibration and the commissioning, biennial calibration has been set at the time of initial delivery. However, however, you will automatically be made aware of the expiry of the calibration interval by ZEISS within one year of commissioning.

- Reset the calibration interval to "annual" after expiry of the first calibration interval.

Calibrating the IC ionization chamber

The use of the ionization chamber IC on the UNIDOS E dosimeter manufactured by PTW requires regular calibration of the ionization chamber IC. PTW recommends calibration every two years.

WARNING:

If the operator does not observe the calibration intervals specified for the IC ionization chamber, the measurements performed during the system check can be incorrect. This can lead to over or under-irradiation.

- Have your IC ionization chamber calibrated every two years in accordance with the recommendations of PTW.

If shorter calibration intervals are stipulated by national codes and regulations, you must comply with them.

Calibration

Initial calibration is performed by ZEISS in the manufacturing plant. As standard practice, only calibrated INTRABEAM Core Systems or single components and calibrated accessories are delivered.

You can perform further calibrations using the application software.

ZEISS Service offers calibration services in the manufacturing plant.


- Comply with the calibration intervals of the XRS X-ray radiation source and the IC ionization chamber, see Page 209.

Activating the function

- Open the "Calibration" screen mask by clicking the <Calibration> tab.

"Calibration" screen mask

You can activate the following functions in the Screen Mask menu "Calibration":

- 1 Date/Time:**
Synchronize the date and time settings of the control unit to those of the user terminal, see Page 212.
- 2 Deflection:**
Ascertain the sensitivity of the beam deflection of a XRS x-ray probe and enter it in the corresponding Calibration file, see Page 214.
- 3 IRM linearity:**
The linear relationships between the individual beam flows and the corresponding IRM count rate is specified at a ratio to the count rate for the maximum beam flow, see Page 218.
- 4 Barometer:**
Compare and correct the measured air pressure in your INTRABEAM Core System with a calibrated, external barometer, and correct, if applicable, see Page 220.

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | Review Treatment | Calibration

5 Thermometer:

Compare and correct the measured temperature of the PAICH in your INTRABEAM Core System using a calibrated, external thermometer, if applicable, see Page 223.

6 Only for old systems: Dilator alignment

Check the straightness of the dilators and center and straighten the dilators, if necessary, see Page 226.

INTRABEAM SYSTEM PRS 500

File Reports Tools About

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | Review Treatment | Calibration

Choose XRS: 000170 40.0 kV

IRM Rate [Hz]: 0

Interlock Status: PDA

Date/Time: 00:00

Deflection	
IRM Linearity	
Barometer	
Thermometer	
Dilator Adjustment	

X beam position [V]: 0.00
Y beam position [V]: 0.00
Target Size [mm]: 2.20

Deflection Calibration

X Static Offset: N/A Volts N/A mm

Y Static Offset: N/A Volts N/A mm

Deflection Factor: N/A Volts/mm

STEP # 0 of 6

Press 'Start' to begin Deflection and Static Offsets Calibration or 'Stop' to terminate.

Date and time synchronization

Activating the function

- Open the "Calibration" screen mask by clicking the <Calibration> tab.

You can synchronize the date and clock settings of the user terminal and PRS500 control console. The PRS500 control console adopts the date and time settings of the user terminal.

The "Date/Time" button (1) is located above the menu in the "Calibration" screen mask.



- Press the "Date/Time" button (1) to activate the "Time & Date" window

which shows the current date and time settings of the user terminal (workstation) and PRS500 control console.

- Press the "Sync" button (2) to start the date and time synchronization.
- Wait until the date and time of the PRS500 control console have been set to the same values as in the user terminal (workstation).

After successful date and time synchronization, the labeling of button (3) changes from <Cancel> to <Done>.

- Press <Done> button (3) to end the routine and close the "Time & Date" window.


(Blank page, for your notes . . .)

Determining the deflection sensitivity

Activating the function

- Open the "Calibration" screen mask by clicking the <Calibration> tab.

You can determine the sensitivity of the deflecting coils of an XRS probe, i.e. you can create a characteristic which shows the voltage change that produces a specific change in the direction of the electron beam in the probe tip. The setup of components and the procedure are the same as for the "Dynamic Offsets" adjustment, see Page 142.

You can then transfer the newly determined deflection coil sensitivity to the calibration file of the relevant XRS X-ray source using the "Deflection Sensitivity History" window. We recommend that you perform this calibration at least once a year and when you receive a new calibration file.

- Press the  "Deflection" button to activate the window required.

Video instruction

Press the  button to start an SQA or system calibration procedure, and the application software will visualize the steps to be performed in a video instruction with displayed text, see Page 136. The application software checks internally whether all verification devices required have been connected and are in position. If this is not the case, the procedure cannot be continued.

- Press the  "Start" button to start an automatic measuring process which will take a few minutes.

After completion of this process, the measured values are displayed in the "Deflection Calibration" window (see the screen mask on the next page), and the "Deflection Factor History" window opens automatically.

You can stop the automatic measuring process by pressing the  "Stop" button.

When the sensitivity of the deflecting coils has been redetermined, you should readjust the beam deflection using the SQA routine "Dynamic Offsets".

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | Review Treatment | Calibration

"Deflection Calibration" screen mask

Consistent system

File Reports Tools About

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | Review Treatment | Calibration

Choose XRS: 000170 40.0 kV

IRM Rate [Hz]: 0

Interlock Status: PDA

Date/Time:

Deflection	
IRM Linearity	
Barometer	
Thermometer	
Dilator Adjustment	

X beam position [V]: 2.21
Y beam position [V]: 0.39
Target Size [mm]: 2.20

Y [mm]

X [mm]

Deflection Calibration

Z PDA Signal [V]

Applied Voltage [V]

X Static Offset: 0.493 Volts 0.312 mm
Y Static Offset: 0.391 Volts 0.247 mm
Deflection Factor: 1.582 Volts/mm

STEP # 6 of 6


Deflection calibration procedure has been completed.

Stop

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | Review Treatment | Calibration

**"Deflection Factor History" screen mask**


1

(Blank page, for your notes . . .)

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. A vertical red margin line runs down the right side of the page. The paper appears to be from a notebook or a standard ruled sheet of paper. There are no markings, text, or drawings on the page.

Determining the linearity of the IRM internal radiation monitor

Activating the function

- Open the "Calibration" screen mask by clicking the <Calibration> tab.
- Press the  "IRM Linearity" button to activate this function.

IRM linearity:

For the IRM, the linear correlations between the individual beam currents with the associated count rates and the count rate at maximum beam current are determined. That means that the linear dependent relationship of the IRM count rate can be determined based on the maximum current (40 μA). An exactly linear behavior would then lead to a count rate of precisely 25% at 10 μA .

Where the 40 kV accelerating voltage is concerned, this step is not applicable because the treatment at 40 kV is only possible with 40 μA .

We recommend that you perform this calibration at least once a year and when you receive a new calibration file.

Only currents producing a minimal count rate are entered in the calibration file and can be used for treatment planning.

Video instruction

Press the  button to start an SQA or system calibration procedure, and the application software will visualize the steps to be performed in a video instruction with displayed text, see Page 136. The application software checks internally whether all verification devices required have been connected and are in position. If this is not the case, the procedure cannot be continued.

- Press the  "Start" button to start an automatic measuring process which will take a few minutes.

After completion of this process, the measured values are displayed in the "IRM Linearity" window (see the screen mask on the next page) and are automatically transferred to the calibration file of the relevant XRS X-ray source.

You can stop the automatic measuring process by pressing the  "Stop" button.

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | Review Treatment | Calibration

**"IRM Linearity" screen mask**

Cross INTRABEAM System

File Reports Tools About

Equipment Inventory | System Quality Assurance | Enter Prescription | Verify Treatment Parameters | Treatment | Review Treatment | Calibration

Choose XRS: 000170 40.0 kV

IRM Rate [Hz]: 0

Interlock Status: PDA

Date/Time:

Deflection	
IRM Linearity	
Barometer	
Thermometer	
Dilator Adjustment	

X beam position [mm]: 0.03
Y beam position [mm]: -0.09
Radial Position Error [mm]: 0.10

IRM rate vs. current

Current, μ A	IRM Count Rate [Hz]	IRM Count Rate / IRM Count Rate @ max μ A	Previous value
40.000	13232.4	1.000	1.000
20.000	6732.0	0.509	0.508
10.000	3361.0	0.254	0.256
5.000	N/A	N/A	***

STEP # 4 of 4

0 25 50 75 100

IRM Linearity Calibration has been completed.

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The INTRABEAM Core System requires the precise air pressure to calculate the correction factor of the dosage during the "PAICH XRS Output Check (Dose rate)", see field in screen mask "Check XRS dose rate", Page 157.

An air pressure sensor is integrated in the control unit. To calibrate this air pressure sensor, you require a calibrated, external barometer. During the calibration, the control unit and the calibrated, external barometer must be located in the same room. Determine the calibration interval according to local legislation.

- Calibrate your air pressure sensor at least once a year.

Activating the function

- Open the "Calibration" screen mask by clicking on the <Calibration> tab.
- Press the "Barometer" button to activate the "Barometer Calibration" window.

The field "Internal Reading" indicates the measured air pressure value for the air pressure sensor integrated in the control unit. As the user of the INTRABEAM Core System, you can carry out a 1 point calibration, see pull-down menu (1).

On activation of the screen mask, <Enter> button (3) is disabled and button (4) has a <Cancel> function. The <Cancel> button permits you to exit the mask without initiating a computation.


- Read off the current air pressure on the calibrated, external barometer and enter it in the field "External reading".
- Use the same measuring unit in the application software and the external barometer (box (2)).

WARNING

Differently set units of measure [mm Hg], [kPa] or [bar] for the air pressure values in the application software and on the calibrated barometer will result in an incorrectly calculated dose rate. This may lead to over-irradiation.

- When entering the air pressure values, please ensure that the same units of measure have been set in the application software and for the calibrated barometer.

When an entry has been made in the "External Reading" box, the <Enter> button (3) is enabled and button (4) has a <Done> function.

- Press the <Enter> button so that the application software computes and enters the offset value for the air pressure sensor integrated in the PRS500 control console.
- Press the <Done> button to close the "Barometer Calibration" screen mask and to transfer the offset value to the INTRABEAM Core System.

Choose XRS: NONE | IRM Rate [Hz]: 0 | Interlock Status: TREATMENT | Date/Time: 14.01.2014 14:00

Barometer Calibration

1 Point Calibration

Select '1 Point Calibration' for periodic calibration of the barometric pressure sensor in the Control Console. Enter the barometric pressure as indicated by a traceable calibrated barometer. Click the 'Enter' button to send the barometric offset correction to the control console.

Internal Reading: 95.7 kPa

External Reading: 95.7 kPa

Enter Done

2 3

$P[\text{kPa}] = \text{Gain}[\text{kPa/V}] * \text{Voltage}[\text{V}] + \text{Offset}[\text{kPa}]$

	Gain[kPa/V]	Offset [kPa]
Values before calibration	110.07	11.89
Values after calibration	110.07	11.83

X beam position [V]: 0.00
Y beam position [V]: 0.00
Target Size [mm]: 0.00

Y [V]
X [V]

"Barometer calibration" screen mask

The pulldown menu (1) offers you the following menu items for selection:

- 1 point calibration,
Can be performed with the procedure described here.
- 2 point calibration,
This function is only for ZEISS Service!
- Load factory settings
Permits you to set your INTRABEAM Core System to a default air pressure value.
- Load last saved value
Permits you to reset your INTRABEAM Core System to the last saved air pressure value.

Barometer Calibration

1 Point Calibration

Select '1 Point Calibration' for periodic calibration of the barometric pressure sensor in the Control Console. Enter the barometric pressure as indicated by a traceable calibrated barometer. Click the 'Enter' button to send the barometric offset correction to the control console.

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
The exact temperature is required by the INTRABEAM core system to calculate the correction factor for the dose during the "PAICH Output Check (Dose rate)", see display field in screen mask "Check XRS Dose Rate", Page 157.

A temperature sensor is integrated in the adapter for the IC Ionization Chamber in the PAICH probe centering unit / IC e Ionization chamber holder. To calibrate this temperature sensor, you require a calibrated, external thermometer.

- Position the temperature probe of the calibrated, external thermometer in the PAICH adapter for the IC ionization chamber. Wait until the temperature reading which the calibrated, external thermometer displays, has stabilized; i.e. wait until the temperature reading hardly changes at all any more.
- Calibrate the integrated temperature sensor at least once a year.

Determination of an incorrect temperature leads to incorrect calculation of the dose value.

- Follow the instructions for calibration of the temperature sensor.


Determine the calibration interval according to local legislation.

Activating the function

- Open the "Calibration" screen mask by clicking on the <Calibration> tab.
- Press the  "Thermometer" button to activate the "Thermometer Calibration" screen mask.

Video instruction

Press the  button to start an SQA or system calibration procedure, and the application software will visualize the steps to be performed in a video instruction with displayed text, see Page 136. The application software checks internally whether all verification devices required have been connected and are in position. If this is not the case, the procedure cannot be continued.

The "Internal Reading" box indicates the reading of the temperature sensor integrated in the PAICH. You can perform a 1-point calibration, see pull-down menu (1).

On activation of the screen mask, <Enter> button (3) is disabled and button (4) has a <Cancel> function. The <Cancel> button permits you to exit the mask without initiating a computation.

- Read off the current temperature for the calibrated, external thermometer. Enter the current temperature reading in the "External Reading" field.
- Use the same measuring unit in the application software and the external thermometer (box (2)).

**WARNING**

Differently set units of measure [Deg C] or [F] for the for the temperature readings in the application software and on the calibrated thermometer will result in an incorrectly calculated dose rate. This may lead to over-irradiation.

- When entering the temperature readings, please ensure that the same units of measure have been set in the application software and for the calibrated thermometer.

If a value has been entered in the "External reading" field, the Enter key (3) becomes active and key (4) has the function <Fertig>.

- Press the <Enter> button so that the application software computes and enters the offset value for the temperature sensor integrated in the PAICH.
- Press <Fertig> to close the "Thermometer calibration" screen mask and to transfer the offset value into the INTRANEAM core system.

"Thermometer Calibration" screen mask

The drop-down menu (1) provides the following options:

- 1 Point Calibration,
can be performed by the user with the procedure described here.
- Load Default
Permits you to set your INTRABEAM Core System to a default temperature value.
- Load Last Saved
Permits you to reset your INTRABEAM Core System to the last saved temperature value.

Choose XRS: NONE

Interlock Status: TREATMENT

IRM Rate [Hz]: 0

Date/Time: [blank]

Deflection: [blank]

IRM Linearity: [blank]

Barometer: [blank]

Thermometer: [blank]

Dilator Adjustment: [blank]

X beam position [V]: 0.00

Y beam position [V]: 0.00

Target Size [mm]: 0.00

X [V]: -4.0 to 4.0

Y [V]: -4.0 to 4.0

Thermometer Calibration

1 Point Calibration

Place a calibrated thermometer into the PAICH in place of the Ionization Chamber. Let the temperature come to equilibrium. Enter the temperature indicated by the calibrated thermometer and click the 'Enter' button.

Internal Reading: 23.50 Deg C

External Reading: 23.50 Deg C

Enter Done

Temperature[DegC] = Gain[DegC/V] * Voltage[V] + Offset[DegC]

	Gain[DegC/V]	Offset[DegC]
Values before calibration	65.02	-29.14
Values after calibration	65.02	-29.82



For predecessor systems only: Dilator Adjustment

Dilators are not part of the INTRABEAM system with PRS 500. The "Dilator Adjustment" function is only described here for customers who have predecessor systems.

Activating the function

- Open the "Calibration" screen mask by clicking the <Calibration> tab.

The setup of components and the procedure for Dilator Adjustment are similar to those used for the "Probe Adjuster" procedure, see Page 138.

- Press the  "Dilator Adjustment" button to activate the "Dilator Adjuster" screen mask.

Troubleshooting table

The INTRABEAM system components are high-quality products. Nevertheless, malfunctions and error messages may occur. The following table, "Help in the event of malfunctions" will assist you if a fault has occurred or an error message appears.

WARNING

The INTRABEAM system does not contain any components that can be serviced by the user. If an unauthorized person tries to service the INTRABEAM system components, uncontrolled X-radiation may be emitted and components damaged. This may lead to over-irradiation.

- If a malfunction occurs or an error message appears which you are unable to rectify using the "Help in the event of malfunctions" table, please label the unit as non-functional and notify ZEISS Service.
- Do not open any components of the INTRABEAM system.
- Only have ZEISS Service or a specialist appointed and authorized by ZEISS perform the service work on your INTRABEAM system.

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Your system has been included in the Carl Zeiss "Service Identification Program" (SIP). A SIP label with a SIP number has been attached to your system, see Page 46.

The SIP number is a unique identification number assigned to your device or system. A product file is maintained by ZEISS for this SIP number.

- For any service requirements, please contact your local contact at the ZEISS subsidiary in your country. You will find their phone number, FAX number and email address on the SIP label.

Fault or Error message	Possible cause	Remedy	See
Communication between control unit and user terminal interrupted.	One or several cables have been detached.	Connect the cables properly.	Page 92
	User terminal (PC) not correctly booted.	Restart the user terminal.	Page 96
Application software no longer runs		- Restart the application software. - If this occurs again, please call ZEISS Service	Page 96
Keyboard does not work	The Y cable or the USB cable is loose or plugged in incorrectly	- Shut down the user terminal. - Connect the Y cable or USB cable properly. - Restart the user terminal.	Page 101 Page 92
IP 10015 You must work OFFLINE because the responding control unit is not designated in the QA Instrument Inventory table. The control unit's response is = xxxx %s. You can log the control unit into the Device Manager, and then log in again in the Main Menu -> Devices -> Login.	The application software has detected that the serial no. of the control unit currently connected does not correspond to the control unit enabled in the Equipment Inventory.	Start with "Offline", in order to set the currently used control unit to "active" in the Components manager.	Page 98
		If the current control unit is not yet included in the Components Manager, then add it and then restart the application software.	Page 112
		Replace the control unit in line with the serial number displayed and connect the RS232 cable.	Page 92
		- Switch on the new control unit and wait until the self-test has been completed successfully. - Click on "Retry" in the Login screen.	Page 96 --

Fault or Error message	Possible cause	Remedy	See
IP 10016 Dosimeter not found! Please ensure that the dosimeter is ON and the serial cable is connected. Select RETRY when ready, or select Continue w/o Dosimeter, in which case the Dosimeter may be turned on later in the session.	Dosimeter switched off or RS232 cable for user terminal is not connected.	• Check the RS232 cable between user terminal and dosimeter and switch on the dosimeter. • Click on "Retry" in the Login screen.	Page 92 --
IP 10018 The responding dosimeter is not "Active" in the QA Instrument Inventory table. The dosimeter response is = xxxx %s. Continue without the dosimeter and enter the dosimeter in the Device Manager, then log in again in the Main Menu -> Options -> Login.	The application software has detected that the serial no. of the dosimeter currently connected does not correspond to the active dosimeter in the Equipment Inventory.	• Start "Offline" and then set the currently used dosimeter to active in the Components Manager. • If the current dosimeter has not yet been entered in the Equipment Inventory, do this now and restart the application software. • Change the dosimeter in accordance with the displayed serial number and connect the RS232 cable. • Switch on the dosimeter. • Click on "Retry" in the Login screen.	Page 92 / Page 96 / Page 112 Page 112 Page 92 Page 96 --
IP 001 Invalid user name or password. Please try again or close the login window.	An invalid password and/or user name was entered.	• Repeat login with correct user name and password.	Page 97
IP 201 Please confirm the calibration data for the XRS x-ray source.	New XRS calibration data has been loaded.	• Click <OK> and check if the serial number and data are correct. • Click <Accept> to enter the new data. • Click <Discard>, to discard the new data.	-- -- --

Fault or Error message	Possible cause	Remedy	See
IP 202 Please confirm the applicator calibration data.	New applicator calibration data has been loaded.	- Click <OK> and check the serial number and data. - Click <Accept> to enter the new data. - Click <Discard>, to discard the new data.	-- -- --
IP 203 Previously existing *.cal file is missing.	Original *.cal file is no longer in the directory originally intended for it.	- Click <OK>. - The corresponding file may be in the Archive folder as a newer *.cal has been loaded. - Or the file has been re-named or deleted.	--
IP 10019 Filament not ready.	XRS has not reached its operating temperature.	Wait until the XRS x-ray source has reached operating temperature. Re-start the procedure, if required.	--
IP 10021 Control console is performing a treatment. Please wait until the treatment is completed. Quit the program or go offline.	User terminal is attempting to communicate with a control unit that is in treatment mode. This is a constellation which has not been anticipated.	- Wait until the treatment session has been completed, and do not use the user terminal. - You can recognize the end of the treatment on the control unit's display or when the "X-Rays on" LED" goes off. - Restart the application software.	Page 191 -- --
IRM Rate error 0x00000008	- System does not reach the expected IRM count rate. - The IRM count rate is reduced due to slight bending of the XRS probe tip.	- Check the serial number of the connected XRS and the XRS selected in the application software. - Based on the diagram of the automatic electron ray centering, check whether the end of the XRS x-ray source is de-centered.	Page 114 Page 218

Fault or Error message	Possible cause	Remedy	See
"0x00000400 High voltage error" "0x00000800 Beam current error" (The error may typically occur briefly after the startup.)	- Internal kV and μ A moni- toring has identified a de- viation.	• XRS X-ray radiation source has not yet as- sumed the correct tem- perature. Allow to achieve the cor- rect temperature. • XRS x-ray radiation source cable not completely plugged in. Insert the ca- ble. • Repeat the procedure.	--
"0x00100000 Control error"	Irradiation has been aborted because the automatic elec- tron beam centering was no longer able to compensate for the decentering of the XRS probe tip.	• Check and rectify the po- sition of the XRS x-ray ra- diation source for mechanical strain and rec- tify if applicable. • If you want to continue the radiation, use the data on the screen mask "Re- view" to plan a new radi- ation treatment for the incorrect dose. In this case, the XRS X-ray radiation source must be re-tested! A different, previously prepared XRS X-ray radi- ation source can also be used.	----- Page 199 Page 130 ff
0x00001000 IRM PIC error (IRM PIC: Monitoring the IRM rate on the main board of the PRS control unit 500). Ex- pected IRM rate outside toler- ance.	Possible defect in the XRS ca- ble connectors or error in the XRS cable or hardware error in control unit or XRS X-ray radiation source.	• Check both cable plug-in connections of XRS cable. • Carry out treatment by choosing function: "Con- tinue same patient".	Page 93 Page 166
0x00004000 VIT PIC error Control console problem during evaluation of the sig- nals	Hardware error in control unit.	• Replace control unit or XRS cable or XRS x-ray ra- diation source.	Page 93

Fault or Error message	Possible cause	Remedy	See
For predecessor systems only:			
0x00002000 ERM PIC error Expected ERM rate exceeds tolerance	Possible faults in cable plug-in connections of ERM cable or fault in ERM cable or hard- ware error in control unit,	• Replace ERM or ERM ca- ble or possibly control unit.	Page 93

Changing the line fuse


The fuse is located in the power switch and power inlet module (1) at the back of the control console.

For the ordering data of replacement fuses, see Page 242.

If required, have the fuse changed by our service staff or an authorized person.

To change the fuse, proceed as follows:

- Turn off the control console using the power switch.
- Disconnect the power cord from the control console.
- Using a flat-tip screwdriver, pry cover (1) up.
- Slide out the fuse holder.
- Remove the defective fuse.
- Insert a new fuse in the fuse holder and slide the fuse holder into the fuse module.
- Close the fuse module with cover (1).
- Connect the power cord to the power socket and switch on the control console using the power switch.

Resterilization

- To protect patients and operating personnel against infections, only use components during the treatment that have been carefully cleaned, disinfected and sterilized.
- Use pre-sterilized, disposable drapes.

WARNING

Sterility cannot be guaranteed or damage can occur if non-permissible methods for cleaning, disinfecting and sterilizing components are used or if components are cleaned, disinfected and sterilized that are not intended for this. This may lead to over-irradiation.

- Only use the permissible cleaning, disinfection and sterilization methods described here.
- Do not clean, disinfect and sterilize any disposable products.

INTRABEAM Components and Accessories	Washing	Rinsing	Disin- fection	Sterile single-use product or recondition- able product
NC32 INTRABEAM Floor Stand	X	-	X	-
INTRABEAM cart	X	-	X	-
PRS 500 control unit	X	-	X	-
User terminal	X	-	X	-
Keyboard	X	X	X	-
XRS 4 X-ray source	X	-	X	-
PDA photo diode array	X	-	X	-
PAICH	X	-	X	-
UNIDOS E dosimeter*	X	-	X	-
Ionization chamber IC*	X	-	X	-
XRS, PDA/PAICH cables	X	X	X	-
V-block	X	X	X	-
X-block	X	X	X	-
INTRABEAM Water Phantom	X	-	X	-
Transport and storage container	X	X	X	-

INTRABEAM Components and Accessories	Washing	Rinsing	Disinfection	Sterile single-use product or reconditionable product
INTRABEAM shields	-	-	-	Sterile Single-use product
INTRABEAM spherical applicators	X	X	X	steam-sterilizable
INTRABEAM Surface Applicator set	X	X	X	steam-sterilizable
INTRABEAM Flat applicator set	X	X	X	steam-sterilizable
INTRABEAM needle applicator	-	-	-	Sterile Single-use product
INTRABEAM drapes	-	-	-	Sterile Single-use product

* If you have to clean and disinfect the UNIDOS E Dosimeter and the IC Ionization Chamber, please refer to the corresponding Instructions for Use by PTW, Freiburg.

Cleaning and disinfecting components

NOTE

If a user with a nickel allergy comes in contact with the XRS housing or a verification device, allergic reactions may occur.

- If required, wear protective clothing and gloves.

CAUTION

The probe tip of the XRS X-ray source contains beryllium (Be). The exterior surface of the probe tip is provided with a protective coating. If the coating is damaged, skin and clothing may be contaminated.

- Avoid any direct skin contact with beryllium and immediately return the XRS X-ray source to the manufacturer.
- First Aid Measures: Wash with lots of water after skin contact with beryllium. Remove contaminated clothing.

NOTE

If the components are not stored in the respective transport or storage case, they can be damaged and contaminated and therefore become unusable.

- Always store the components in their transport and sterilization case when you are not using them.

NOTE

If the interlock surfaces are contaminated or have become matte, the interlock may not be satisfied and the INTRABEAM system may not generate X-radiation.

- Use only approved cleaning agents to clean the interlock surfaces.

The procedure described below applies to the following components of the INTRABEAM Core System:

- XRS X-ray source and probe guard
- PRS 500 control unit
- Cables (XRS, PAICH/PDA connecting cables)
- V-block guide and X-block
- Transport and storage case for PDA, PAICH, and XRS
- PAICH probe adjuster/ionization chamber holder
- PDA Photo Diode Array
- Please perform the following steps:
 - Thoroughly clean all components as soon as possible after use and prior to disinfection.
 - Check all surfaces of each item for gross contamination and impurities.
 - Remove contamination or impurities using a single-use sponge, or a gauze pad moistened with water.
 - Continue washing the surfaces until all visible contaminants are removed. If necessary, a non-abrasive cleaning cloth may be used.
 - Disinfect each item using a single-use sponge or gauze pad moistened with a solution of one part bleach (4-7% hypochloride) to ten parts distilled water (by volume).
 - For difficult-to-reach areas (such as the inside of tubes), substitute a lint-free swab or soft-bristled brush for the gauze pad.
 - Allow the items to dry.

Cleaning the touchscreen

Deactivate the touch function to avoid unintended faulty operation while cleaning the touchscreen, see Page 62:

- Press the <Deactivate touch function> button.
When the touch function has been deactivated, the "Touch" indicator lamp illuminates.
Press this button again to reactivate the touch function.

NOTE

The wrong cleaning method can damage the user terminal, e.g. the paint can peel off. The wrong cleaning method can also lead to malfunctions and unexpected behavior of the user terminal.

- Use only the recommended cleansers and cleaning methods for the user terminal:
 - Do not spray any cleanser or disinfectant onto the user terminal or system carrier as moisture may seep into the devices.
 - Do not use solvents or caustic and abrasive cleansers or disinfectants to clean the mechanical surfaces.
- Clean the touchscreen with alcohol-free glass cleaner.
- Apply the cleanser to a soft, clean cloth (do not spray or pour it directly on the touchscreen!) and wipe the touchscreen with the moist cloth.

Cleaning the equipment trolley(INTRABEAM Cart)

- If possible, clean the INTRABEAM cart and the instruments immediately after use. Contaminations should not be allowed to dry on the objects, as this would make cleaning and disinfecting more difficult.
- Wipe the INTRABEAM cart with a moist cloth. Remove any residue using a mixture of 50% ethyl alcohol and 50% distilled water plus a dash of household dish-washing liquid.

NOTE

The use of unsuitable cleaning and disinfection agents may lead to damage and corrosion of the INTRABEAM cart surface.

- Do not use any aggressive or abrasive cleaning agents.

Disinfecting user terminal and cart

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The maximum concentrations are:

- For alcohol (tested with 2 propanol): 60%
- For aldehyde (tested with glutaraldehyde): 2%
- For quaternary compounds (tested with DDAC): 0.2%

Casters of the system cart

- Clean and disinfect the casters of the system cart using a standard disinfectant.

NOTE

The casters of the system cart will be damaged if they are left for too long in a disinfection bath.

- Give the disinfectant time to take effect and then wipe it off the casters.

Handling, transport and shipment

- Wear protective clothing such as gloves, an apron or gown and an eye shield when handling contaminated components of the INTRABEAMCore system.
- Always assume INTRABEAM Core System components to be contaminated unless they are known to be clean, disinfected and sterile.
- When transporting the INTRABEAM Core System or any of its components, take appropriate care to protect the items being transported against shock and severe variations in temperature.
- When transporting INTRABEAM Core System components inside the hospital, always use the transport and sterilization containers provided for this. Exceptions are the control unit, the user terminal and the UNIDOS dosimeter.
- If transporting INTRABEAM Core System components outside the hospital, they must always be in a clean and disinfected condition and be in the original packaging or the shuttle containers provided by ZEISS.
- Always ship INTRABEAM Core System components in a clean and disinfected condition using the original packaging or the shuttle containers provided by ZEISS.

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- In order to dispose of outer packaging materials please refer to the local guidelines and bylaws.

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(One of each for the NC32 INTRABEAM Floor Stand and the isolating transformers in the INTRABEAM cart)

Europe	000000-0135-487
USA	397902-9066-000
UK	309850-9003-000
Switzerland	000000-0426-193
Argentina	000000-0434-527
China	000000-0475-509
Brazil	000000-0527-729

Spare parts

- Only use the spare parts included in the delivery package.
- Use the spare parts included in the delivery package only for the intended purpose described.
- Please send your orders to your local contact of the ZEISS subsidiary in your country.

Designation	Order number
Line fuse T1.25A E/250V; 5 x 20 mm for PRS500 control console, for (100-240) VAC, 1 pce	000000-0302-911 (2 pcs required)
UNIDOS E Dosimeter (complete unit)	000000-0485-223
Line fuse 315 mA T (1 pc), for UNIDOS E dosimeter, for (100 - 240 VAC	000000-0127-015 (2 pcs required)
Interlock control plug	000000-1890-246

Accessories

- Only use the accessories included in the delivery package.
- Use the accessories included in the delivery package only for the intended purpose described.
- Please send your orders to your local contact of the ZEISS subsidiary in your country.

Designation	Order number
Cable (to RS 232)	000000-0466-623
CD-R blanks (pack of 50)	000000-1195-144

Spherical Applicators

You can order the spherical applicators either as single items or as a set. The set comprises eight applicators with different sphere diameters and a sterilization container. The spherical applicators and the sterilization container are also available as single items.

Designation	Order number
INTRABEAM spherical applicator set (8 applicators in a sterile container)	304534-6000-348
INTRABEAM spherical applicator Ø1.5 cm	304534-6000-557
INTRABEAM spherical applicator Ø2 cm	304534-6000-558
INTRABEAM spherical applicator Ø2.5 cm	304534-6000-559
INTRABEAM spherical applicator Ø3 cm	304534-6000-560
INTRABEAM spherical applicator Ø3.5 cm	304534-6000-561
INTRABEAM spherical applicator Ø4 cm	304534-6000-562
INTRABEAM spherical applicator Ø4.5 cm	304534-6000-563
INTRABEAM spherical applicator Ø5 cm	304534-6000-564
Sterile container with reusable filter (also used as transport and storage case)	304534-6000-347
Reusable filter	304534-3000-167

Radiation shields

Designation	Cat. No.
Radiation shield, flat	304534-6000-595

Drape

INTRABEAM drape (sterile covers for the NC32 INTRABEAM Floor Stand)	326090-0000-000
--	-----------------

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System data

Technical data of the XRS 4 x-ray source	246
Technical data for UNIDOS E	247
Technical data of the PRS 500 control console.....	248
Technical data of the user terminal.....	249
 External interlock and X-radiation indicator	250
External X-ray interlock.....	250
External X-ray display	251
 EMC (electromagnetic compatibility)	252
Electromagnetic radiation disturbance	254
Electromagnetic immunity	255
Recommended safety distances	258
 Environmental conditions	259
 Regulatory information	260
Regulatory information on the PRS 500 controller	260
Regulatory Information about the User Terminal	261

Technical data of the XRS 4 x-ray source

Mechanics

Weight	1.62 kg
Dimensions	(175 x 110 x 70) mm
Probe	3.2 mm diameter, length 100 mm coated with chromium nitride (CrN)

Radiation output

Dosage characteristics

All stated values are valid for operation at maximum settings of 50 kV and 40 μ A.

10	9.8	1.5
20	2.3	6.5
30	0.9	16.7
40	0.4	37.5
50	0.25	60

Radiation output stability

$\pm 5\%$ for first 5 minutes of X-ray output; then $\pm 2\%$.

Radiation distribution

- Point-source type emission
- Dose isotropic about XRS isocenter
- Sharp dosage drop (approx $1/r^3$) in soft tissue or equivalent.
- Target accuracy of the delivered dose better than ± 1 mm for treatment diameter of 40 mm.

Self-filtration

50	0.64	20.4
40	0.48	19.1

Radiation exposure

Radiation exposure due to an exposed XRS x-ray source at 50 kV and 40 μ A at a distance of 1 meter equals 9R/hour. (1.32 mGy/min).

Worst case radiation leakage from a properly shielded XRS x-ray source with 50 kV and 40 μ A at 50 mm is 0.6 mR/hour. (0.0035 mGy/min).

Standard IEC 60601-2-8 allows a radiation exposure at a distance of 50 mm for radiation leakage of ≤ 3 mR/hour. (≤ 0.134 mGy/hour).

Technical data for UNIDOS E

Weight	3.2 kg (7.05 lb)	
Dimensions (W x H x D)	250 mm X 100 mm X 260 mm (9.84" X 3.94" X 10.24")	
Start-up time	5 min	
Power rating	Input voltage	(100 - 230) V $\pm$ 15%
	Rated frequency	(50 - 60) Hz
	Current consumption	0.05 A

Technical data of the PRS 500 control console

Weight	4.5 kg (9.92 lb)
Dimensions (W x H x D)	(305 x 89 x 381) mm / (12.00 x 3.50 x 15.00)"
Power connection	Only connect the system to power outlets which are equipped with a properly connected protective ground conductor.
Input voltage	100 V - 240 V
Rated frequency	50 Hz - 60 Hz
Power consumption	60 VA
Selectable operating parameters:	
Accelerating voltage /	40 kV / 40 μ A
Beam current	50 kV/ 5, 10, 20 or 40 μ A
	40 μ A is the maximum X-ray tube current over the entire voltage range.

The system is designed for continuous operation (standby).

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Technical data of the user terminal

PRS 500 user terminal

Weight	12 kg (2.20 lb)
Dimensions (W x H x D)	(400 x 341 x 100) mm / (15.75 x 43.31 x 3.94)".
Power connection	Only connect the device to power outlets which are equipped with a properly connected protective ground conductor.
Input voltage	100 V - 240 V
Rated frequency	50 Hz - 60 Hz
Power consumption	max. 100 VA

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The external interlock connector "EXT INTERLOCK" on the rear panel of the control console has the following two functions:

- Connection facility for a separate switch to be used to stop, prevent or enable the generation of X-ray radiation.
- Connection facility for an additional external indicator (warning lamp) of the generation of X-ray radiation.

The interlock plug attached to the control console matches the "EXT INTERLOCK" connector. Pin 1 and pin 2 of this interlock plug are connected via a wire jumper. When inserted into the connector "EXT INTERLOCK", this interlock plug closes the external interlock circuit.

If you want to use a separate switch, for example, a door switch, this switch must be connected to pins 1 and 2 of connector "EXT INTERLOCK" via a separate switch. This connector must be provided by the customer.

The control console is supplied with a matching plug for the external interlock connector. Plugs matching the external interlock connector can be purchased as spare parts from ZEISS Service, see Page 242 "Interlock Control Plug" spare parts.


Figure to the left: Pin assignment of connector "EXT INTERLOCK".

External X-ray interlock

Place a low impedance (<10 ohm) connection between pins 1 and 2 of the connector.


CAUTION

The system will be damaged if line voltage is conducted by mistake to the PRS 500 control console via the external interlock. The user may come in contact with live line voltage and suffer burns and ventricular fibrillation as a result.

- Do not apply an external voltage to pin 1 or 2 of the connector.

Recommended circuit to be connected to pins 1 and 2:

External X-ray display

The customer is responsible for the supply of the external X-ray radiation display.

ZEISS recommends using an ohmic resistor (R) at pins 3 and 4 to limit the current and a relay to activate an external X-ray radiation display.

Pins 3 and 4 form a closed circuit when X-ray radiation is being generated and an open circuit at all other times.


Maximum voltage on pins 3 and 4

50 V DC or 24 V AC

Maximum current at pins 3 and 4

0.5 A rms

Example:

The connection diagram shown above is only a schematic representation. Connection is only permitted if the requirements of IEC 60601-1-1 2000 or IEC 60601-1:2005, Section 16 (safety of electrical medical systems) are met and the installation has been properly installed. You must also ensure that double insulation in accordance with IEC 60601-1 is provided for the voltage source and for the separation between line power and low voltage sources. National regulations must be observed for the installation.

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Electrical medical devices are subject to special precautionary measures with regard to electromagnetic compatibility (EMC).

The INTRABEAM System PRS 500 with XRS 4 fulfills the EMC requirements of EN 60601-1-2:2007.

- Only install and operate the INTRABEAM system according to the following precautionary measures with regard to EMC.
- Only use components and spare parts approved by ZEISS for this INTRABEAM system.

Additional equipment which you connect to electrical medical devices must demonstrably comply with the applicable IEC or ISO standards (e.g. IEC 60950 for data processing equipment). In addition, all configurations must meet the normative requirements for electrical medical systems (see IEC 60601-1-1:1999 or Section 16 of IEC 60601-1:2005). A person connecting additional equipment to medical electrical devices is the system configurator and is therefore responsible for ensuring that the system complies with the normative requirements for systems.

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- Do not use any portable or mobile RF communication equipment in the vicinity of the INTRABEAM system, as it is not possible to exclude the possibility that the function of the INTRABEAM system will be affected.
- Do not use any cellphones in the vicinity of the INTRABEAM system, as these may represent a potential hazard for the correct function of medical/electrical devices. Malfunctions may occur depending on a variety of local factors. These cannot be foreseen or assessed in any whatsoever.
- Please note the EMC guidelines on the following pages.

The electromagnetic compatibility was ascertained on a fully configured PRS500 INTRABEAM system with XRS 4.

The full configuration of a PRS500 INTRABEAM system with XRS 4 is comprised of the following components:

- INTRABEAM cart REO-MED 600 isolating transformer
Order No.: 000000-1501-288
- INTRABEAM PRS 500 control unit
Order No.: 304534-7000-338

- INTRABEAM PRS 500 user terminal (including keyboard)
Order No.: 304534-7000-352
- NC32 INTRABEAM Floor Stand,
Order No.: 000000-1096-945
- INTRABEAM XRS 4 x-ray source
Order No.: 304534-7000-001
- UNIDOS E dosimeter
Order No.: 000000-0485-223

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- Keyboard with user terminal,
Order No. (e.g. English USB keyboard): 000000-0545-564
- PRS 500 user terminal including PRS 500 control unit:
RS232 connection cable, length 2.5 m
Order No.: 000000-0466-623
- PRS 500 user terminal with UNIDOS E dosimeter:
RS232 connection cable, length 2.5 m
Order No.: 000000-0466-623
- PRS 500 control unit including NC32 INTRABEAM Floor Stand,
i.e. control cable for XRS 4 x-ray source, length 9 m
Order No.: 304534-7000-092
- PRS 500 control unit including PDA/PAICH: RS232 connection cable,
length 1.5 m, Order No.: 304534-7000-091
- Power supply for INTRABEAM Cart, i.e. power connection for
REO-MED 600 isolating transformer: RS232 connection cable, length 6 m
- Power supply for NC32 INTRABEAM floor stand:
RS232 connection cable, length 6 m

**WARNING:**

If you illicitly use components which are not approved by ZEISS for use with the INTRABEAM System (i.e. ME devices, accessories, isolating transformers (transducers) and cables), this can cause increased emissions of electromagnetic interference or a reduced electromagnetic interference resistance in your INTRABEAM system.

A safe, interference-free operation of your INTRABEAM system is no longer possible which may result in the patient receiving excessive levels of radiation.

- Only operate components approved by ZEISS in your INTRABEAM system.

Electromagnetic radiation disturbance

The INTRABEAM system is intended for operation in an electromagnetic environment as specified below. The customer or the user of the INTRABEAM system is responsible for ensuring that the system is operated in such an environment.

Radiation disturbance measurements	Compliance	Electromagnetic environment - guidelines
RF emissions as per CISPR 11	Group 1	The INTRABEAM system uses RF energy only for its internal functions. As a result, RF emissions are very low and unlikely to cause any interference in nearby electronic devices.
RF emissions as per CISPR 11	Class B	The INTRABEAM system is suitable for use in all facilities including locations in residential environments and those directly connected to the public power supply network which also supplies buildings used for residential purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Emission of voltage fluctuations/ flicker as per IEC 61000-3-3	Compliance	

Electromagnetic immunity

The INTRABEAM system is intended for operation in an electromagnetic environment as specified below. The customer or the user of the INTRABEAM system is responsible for ensuring that the system is operated in such an environment.

Immunity tests	Test level as per IEC 60601-1-2	Compliance level	Electromagnetic environment - guide-lines
Electrostatic discharge (ESD) as per IEC 61000-4-2	±2, 4, 6 kV contact discharge ±2, 4, 8 kV air discharge	±2, 4, 6 kV contact discharge ±2, 4, 8 kV air discharge	Floors should be made of wood or concrete or be covered with ceramic tiles. If the flooring contains synthetic materials, the relative humidity must be at least 30%.
Fast transient/burst immunity as per IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	The quality of the supply voltage should be that of a typical business or hospital environment.
Surges as per IEC 61000-4-5	±0.5kV / 1 kV line to line voltage ±0.5, 1, 2 kV line to ground voltage	±0.5kV, 1 kV line to line voltage ±0.5, 1, 2 kV line to ground voltage	The quality of the supply voltage should be that of a typical business or hospital environment.
Short interruptions and voltage variations as per IEC 61000-4-11	5% U_T for 0.5 cycle 40% U_T for 5 cycles 70% U_T for 25 cycles	5% U_T for 0.5 cycle 40% U_T for 5 cycles 70% U_T for 25 cycles	The quality of the supply voltage should be that of a typical business or hospital environment. If the user of the INTRABEAM system requires continued function even in the event of interruptions in the power supply, we recommend to power the INTRABEAM system from an uninterruptible power supply.
Voltage dips as per IEC 61000-4-11	5% U_T for 250 cycles	5% U_T for 250 cycles	
Power frequency (50/60 Hz) magnetic field as per IEC 61000-4-8	50 Hz 3 A/m 60 Hz 3 A/m	50 Hz 3 A/m 60 Hz 3 A/m	Power frequency magnetic fields should be at levels typical of business and hospital environments.

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The INTRABEAM system is intended for operation in an electromagnetic environment as specified below. The customer or the user of the INTRABEAM system is responsible for ensuring that the system is operated in such an environment.

Immunity tests	Test level as per IEC 60601	Compliance level	Electromagnetic environment - guidelines
			Portable and mobile radio communication equipment should not be used closer to the INTRABEAM system, including cables, than the recommended safety distance that is calculated using the equation applicable to the transmission frequency involved.
Conducted RF disturbances as per IEC 61000-4-6	3 V _{effective} value 150 kHz to 80 MHz	[U ₁] = 3 V _{eff} 150 kHz to 80 MHz 80% AM	Recommended safety distance: $d = 1.17 \sqrt{P}$
Radiated RF disturbances as per IEC 61000-4-3	3 V/m 80 - 2500 MHz	[E ₁] = 10 V/m 80 - 2500 MHz 80% AM	Recommended safety distance: $d = 0.35 \sqrt{P}$ for 80 MHz - 800 MHz $d = 0.7 \sqrt{P}$ for 800 MHz - 2.5 GHz

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Field strengths from stationary RF transmitters, as determined by a site survey^a, should be less than the compliance level in all frequency ranges.^b

Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is influenced by absorption and reflection by structures, objects and persons.

^a Field strengths of stationary transmitters such as base stations for mobile telephones and mobile land radio equipment, amateur radio stations, AM and FM radio broadcast and TV broadcast transmitters cannot be theoretically predicted accurately. To assess the electromagnetic environment with respect to stationary RF transmitters, a site study of the electromagnetic phenomena should be considered. If the measured field strength in the location where the INTRABEAM system is used exceeds the compliance levels indicated above, the INTRABEAM system should be monitored to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the ME equipment or ME system.

^b Field strengths should be less than 3 V/m over the frequency range from 150 kHz to 80 MHz.

Recommended safety distances

The INTRABEAM system is intended for use in an electromagnetic environment in which RF disturbances are controlled. The customer or the user of the INTRABEAM system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the INTRABEAM system, depending on the output power of the communication equipment as specified below.

Rated power of the transmitter W	Separation distance depending on transmission frequency m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = \left[\frac{3,5}{U_1} \right] \sqrt{P}$	$d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$	$d = \left[\frac{7}{E_1} \right] \sqrt{P}$
	For $[U_1] = 3 \text{ V}_{\text{eff}}$	For $[E_1] = 10 \text{ V/m}$	For $[E_1] = 10 \text{ V/m}$
0.01	0.12	0.04	0.07
0.1	0.37	0.11	0.22
1	1.17	0.35	0.70
10	3.69	1.11	2.21
100	11.67	3.50	7.00

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NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is influenced by absorption and reflection by structures, objects and persons.

Environmental conditions

For operation	Temperature	+ 15 °C ... + 40 °C + 59 °F ... + 104°F
	Rel. Humidity	30 % ... 75 %
	Air pressure	800 hPa ... 1060 hPa
For transport and storage	Temperature	- 20 °C ... + 70 °C - 4 °F ... + 158 °F
	Rel. Humidity (without condensation)	10 %...90 %
	Air pressure	500 hPa...1060 hPa

Environmental conditions differing from the core system

UNIDOS E, for transport and storage	Temperature	- 20°C ... + 60°C - 4 °F ... 140°F
	Rel. Air humidity (without condensation)	10% ... 85%
	Air pressure	600 hPa ... 1200 hPa

Regulatory information

Regulatory information on the PRS 500 controller

Requirement	Labelling
Electrical safety	The medical device meets the requirements according to – IEC 60601-1:2005, – IEC 601-1:1988+A1:1991+A2:1995, – CAN/CSA-C2.2 No. 601.1-M90. Classification by level of protection against electrical shock: Protection class 1 Protection type (according to IEC 60529): IP 20
Electromagnetic compatibility	The medical device is allocated to class B (acc. to IEC 61000-3-2) in line with IEC 60601-1-2:2007.
CE label	The medical device meets the essential requirements stipulated in Appendix I of the Medical Device Directive 93/42/EEC. In line with Appendix XII of the Medical Device Directive 93/42/EEC, the medical device is labeled with 0297 .

Regulatory Information about the User Terminal

Electrical safety	The medical device meets the requirements according to – IEC 60601-1:2005, – IEC 601-1:1988+A1:1991+A2:1995, – CAN/CSA-C2.2 No. 601.1-M90. Classification by level of protection against electrical shock: Protection class 1 Protection type (according to IEC 60529): IP 20
Electromagnetic compatibility	The medical device is allocated to class B (acc. to IEC 61000-3-2) in line with IEC 60601-1-2:2007.
CE label	The medical device meets the essential requirements stipulated in Appendix I of the Medical Device Directive 93/42/EEC. In line with Appendix XII of the Medical Device Directive 93/42/EEC, the medical device is labeled with 0297 .

1. The first part of the document is a title page. It contains the title of the document, the author's name, and the date of the document. The title is "The History of the United States of America". The author is "John Adams". The date is "1776".

2. The second part of the document is the preface. It contains the author's introduction to the document. The author states that the document is a history of the United States of America, and that it is written for the purpose of informing the public of the events that have shaped the nation.

3. The third part of the document is the main body of the text. It contains the author's account of the events that have shaped the United States of America. The author begins with the founding of the nation, and then proceeds to describe the various events that have shaped the nation's history, including the American Revolution, the Civil War, and the Reconstruction era.

4. The fourth part of the document is the conclusion. It contains the author's final thoughts on the history of the United States of America. The author states that the United States of America is a nation that has been shaped by the events of its history, and that it is a nation that is still being shaped by the events of its present.

5. The fifth part of the document is the index. It contains a list of the events that are mentioned in the document, and the pages on which they are mentioned. This part of the document is useful for finding specific information about the events of the United States of America's history.

Supplementary documents

Applicable user manuals	264
Equipment checklist for the operating room	265
Literature.....	266
Radiotherapy	266

Applicable user manuals

- INTRABEAM cart
Order Number, ZEISS, 000000-1501-288,
ZEISS Instructions for Use G-30-1704
- NC32 INTRABEAM Floor Stand,
Order Number, ZEISS, 1096-945
ZEISS instructions for Use G-30-1476
- INTRABEAM drape
Order Number, ZEISS, 32090-0000-000,
ZEISS instructions for Use G-30-1526
- Radiation Shield, flat
Order Number, ZEISS, 304534-6000-595,
ZEISS instructions for Use G-30-1880
- INTRABEAM spherical applicators
Order Number, ZEISS, 304534-6000-348,
ZEISS instructions for Use G-30-1478
- INTRABEAM Needle Applicator,
Order Number, ZEISS, 304534-3001-200,
ZEISS instructions for Use G-30-1723
- INTRABEAM Surface Applicator Set 1.0 cm, 2.0 cm, 3.0 cm and 4.0 cm
Order Numbers, ZEISS, 304534-3001-530, -531, -532, -533;
ZEISS instructions for Use G-30-1858
- INTRABEAM Flat Applicator Set 1.0 cm, 2.0 cm, 3.0 cm, 4.0 cm,
5.0 cm and 6.0 cm;
Order numbers, ZEISS, 304534-3001-565, -566, -567, -568,
-569 and -570;
ZEISS instructions for Use G-30-1880
- UNIDOS E dosimeter from PTW, Type 10010,
Instructions for Use_en, PTW, 54513100_02
Instrucciones de manejo_es, PTW, 54519600_00
Gebrauchsanweisung_de, PTW, 54515100_02
- INTRABEAM Water Phantom,
Order Number, ZEISS, 304534-7000-600,
ZEISS instructions for Use G-30-1767
- PRS 500 user terminal
Instructions for Use from PENTA
Medical ETX-CD Zeiss manual English R 1.1 (on CD-ROM)

Equipment checklist for the operating room

Item

PRS 500 control console with power cord

User terminal with power cord and keyboard

XRS set:

XRS (X-ray source)

XRS cable

Transport and sterilization tray for XRS set

XRS accessory set:

PAICH (probe adjuster/ionization chamber holder)

PDA (photodiode array)

PAICH/PDA cable

Transport and storage container for XRS accessory set

Ionization Chamber Set

IC (Ionization Chamber)

IC extension cable

IC shield

Dosimeter with power cord

V-block

X-block

Sterile applicators

User Manuals

X-ray survey meter

(not part of the INTRABEAM Core System - recommended)

Rule (not part of the INTRABEAM Core System - recommended)

Literature

A few examples of literature on the INTRABEAM system are listed below. You will find further publications in the list of references on the INTRABEAM system. The list of publications can be ordered from ZEISS.

Radiotherapy

Dosimetry and Physics

- Armoogum KS, Parry JM, Souliman SK, Sutton DG, Mackay CD: Functional intercomparison of intraoperative radiotherapy equipment - Photon Radiosurgery System. *Radiat Oncol* 2007;2:11.
- Armoogum K, Parry J, Mackay C, Thompson A, Vaidya J, Ackland C, Gardner J, Munro A, Eljamel S: Implementation and experiences of an intraoperative radiotherapy service. *Clin Oncol (R Coll Radiol)* 2007;19(3 Suppl):S8-9.
- Ebert MA, Carruthers B: Dosimetric Characteristics of a Low-kV intra-operative X-ray Source: Implications for use in a Clinical Trial for Treatment of Low-risk Breast Cancer. *Med Phys* 2003;30:2424-31.
- Ebert MA, Carruthers B, Lanzon PJ, Haworth A, Clarke J, Caswell NM, Siddiqui SA: Dosimetry of a Low-kV Intra-operative X-ray Source using Basic Analytical Beam Models. *Australas Phys Eng Sci Med* 2002;25:119-23.
- Clarke J, Lanzon P, Haworth A, Ebert M, Caswell N, Tuchnya T: The PeC Intrabeam Miniature X-Ray Device - Physics Implementation. *Proceedings Intercontinental Workshop and Symposium on IORT using Mobile Systems*, Sydney, 2002.
- Biggs PJ, Beatty J, Yasuda T: Distance-dose Curve for a Miniature X-ray Tube for Stereotactic Radiosurgery using an Optimized Aperture with a Parallel-Plate Ionization Chamber. *Med Phys* 1999;26:2550-54.
- Yasuda T, Beatty J, Biggs PJ, Gall K: Two-dimensional Dose Distribution of a Miniature X-ray Device for Stereotactic Radiosurgery. *Med Phys* 1998;25:1212 - 1216.
- Hakim R, Zervas NT, Hakim F, Butler WE, Beatty J, Yanch JC, Biggs PJ, Gall KP, Sliski AP: Initial Characterization of the Dosimetry and Radiobiology of a Device for Administering Interstitial Stereotactic Radiosurgery. *Neurosurgery* 1997;40:510-7.

- Yanch JC, and Harte KJ: Monte Carlo Simulation of a Miniature, Radiosurgery X-ray tube using the IST 3.0 Coupled Electron-Photon Transport Code. *Med Phys* 1996;23:1551-8.
- Dinsmore M, Harte KJ, Sliski AP, Smith DO, Nomikos PM, Dalterio MJ, Boom AJ, Leonard WF, Oettinger PE, Yanch JC: A New Miniature X-ray Source for Interstitial Radiosurgery: Device Description. *Med Phys* 1996;23:45-52.
- Beatty J, Biggs PJ, Gall K, Okunieff P, Pardo FS, Harte KJ, Dalterio MJ, Sliski AP: A New Miniature X-ray Source for Interstitial Radiosurgery: Dosimetry. *Med Phys* 1996;23:53-62.
- Biggs DS and Thomson ES: Radiation properties of a miniature X-ray device for radiosurgery. *Br J Ophthalmol* 1996 69: 544-7

Radiobiology

- Herskind C, Kraus-Tiefenbacher U, Keller A, Fleckenstein K, Wenz F: Intra-operative Boost in mammary cancer using 50 kV X-Rays from INTRA-BEAM: Radiobiological modeling of late reaction. *Rev. Cancer* 2008; 22. No Extra:6.
- Belletti B et al., Targeted Intraoperative Radiotherapy Impairs the Stimulation of Breast Cancer Cell Proliferation and Invasion Caused by Surgical Wounding, *Clin Cancer Res* 1325 2008;14(5)
- Gardner J, Ackland CY, Thompson AM, Vaidya SJ, Dewar JA, Adamson DJ, Armoogum K, Mackay CD, Munro AJ: On the Practicalities of Translational Radiobiology: Lessons from IORT (Intra-operative Radiotherapy). *Clin Oncol (R Coll Radiol)* 2007;19(3 Suppl):S29.
- Herskind C, Schalla S, Hahn EW, Höver KH, Wenz F: Influence of different dose rates on cell recovery and RBE at different spatial positions during protracted conformal radiotherapy. *Radiat Prot Dosimetry* 2006;122:498-505.

Breast Cancer

- Vaidya J et al. Targeted intraoperative radiotherapy (TARGIT) yields very low recurrence rates when given as a boost, *Int. J. Radiation Oncology Biol. Phys.*, Vol. 66, No. 5, pp. 1335-1338, 2006
- Vaidya J et al. Targeted intraoperative radiotherapy (TARGIT) Boost after breast conserving surgery results in a remarkably low recurrence rates in a standard risk population: 5 year results, *Proc. San Antonio Breast Cancer Conference* Dec 2009
- Wenz F et al. Intraoperative radiotherapy (IORT) as a boost during breast conserving surgery (BCS) using low kV X-rays - the first five years of experience with a novel approach, *Int. J. Radiation Oncology Biol. Phys.* (Article 18522)
- Kraus-Tiefenbacher U et al. Postoperative seroma formation after intraoperative radiotherapy using low-kilovoltage X-rays given during breast conserving surgery, *Int. J. Radiation Oncology Biol. Phys.*, 2009 Oct 14.
- -Tiefenbacher U et al. Postoperative seroma formation after intraoperative radiotherapy using low-kilovoltage X-rays given during breast conserving surgery, *Int. J. Radiation Oncology Biol. Phys.*, Vol. 66, No. 2, pp. 377-381, 2006
- Keshtgar M et al. Breast Cosmetic outcome 2 and 3 years after intraoperative radiotherapy compared with external beam radiotherapy for early breast cancer: An objective assessment of patients from a randomized controlled trial, *Breast Cancer symposium* 2009; A 294
- Kraus-Tiefenbacher U et al. Intraoperative Radiotherapy (IORT) as a Boost in Patients with Early-Stage Breast Cancer - Acute Toxicity, *Onkologie*. 2003 Dec;26(6):596-8
- Keshtgar M et al. Intraoperative Radiotherapy for Breast Cancer in patients where external beam radiation was not possible: Results after 3 years of follow up, *Proc. St Gallen Breast Cancer Conference* Mar 2009
- Kraus-Tiefenbacher U et al. Intraoperative radiotherapy (IORT) is an option for patients with localized breast recurrences after previous external-beam radiotherapy, *BMC Cancer* 2007, 7:178

Skin Cancer

- Bodner W et al. Use of low energy X-rays in the treatment of superficial nonmelanomatous skin cancers, *Cancer Investigation* Vol. 21, No. 3, pp. 355-362, 2003

Primary Brain Tumors and Metastases

- Eljamel S et al. Intraoperative Radiotherapy using mobile miniature X-Ray Source in malignant brain tumors, Proc. American Association of Neurosurgery (AANS) 2008
- Kalapurakal JA, Goldmann S, Stellpflug W, Curran J, Sathiaselalan V, Marymont MH, Tomita T: Phase I study of intraoperative radiotherapy with photon radiosurgery system in children with recurrent brain tumours: preliminary report of first dose level (10 Gy). Int J Rad Oncol Biol Phys 2006;65:800-8.
- Curry W et al. Stereotactic interstitial radiosurgery for cerebral metastases, J Neurosurg 103:630-635, 2005

Gastric and Rectal Cancer

- Lyadov K, Sinyakin S, Krymskiy VA, Krymskiy AV: Improvement of curativity of video-assisted surgery for colorectal cancer due to intraoperative contact radiotherapy using the INTRABEAM system. Rev. Cancer 2008; 22. No Extra:8.
- Algur E. et al. Interstitial Photon Radiosurgery System for recurrent and locally advanced rectal cancer: Retrospective review of 24 patients, Proc ASCO GI Symposium 2005
- Lyadov KV, Krymskiy VA, Markin AY, Krymskiy AV: Evaluation of safety of intraoperative radiotherapy using the INTRABEAM system in combined treatment of gastric cancer. Rev. Cancer 2008;22. No Extra:25

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Index

Symbols

"Barometer calibration" screen mask	222
"Calibration" screen mask	210
"Deflection Calibration" screen mask	215
"Deflection Factor History" screen mask	216
"Enter Prescription" screen mask	163
"Equipment Inventory" screen mask	112
"IRM Linearity" screen mask	219
"System Quality Assurance" screen mask	131
"Thermometer Calibration" screen mask	225
"Treatment" screen mask	181
"Verify Treatment Parameters" screen mask	174

A

Abbreviations in these Instructions for Use	3
Accelerating voltage, PRS control unit	248
Acceleration voltage	22
Accessories	51, 243
Accessories for checking the system	14
Accessories for the INTRABEAM core system	80
Accessories of the INTRABEAM system	16
Accessories, shields	244
Additional equipment	24, 252
Address of manufacturer	2
Adjusting (straightening) XRS probe and dilators	70
Administrative tools	103
Administrative tools, user management	97
Air pressure sensor	53
Air pressure sensor, calibration	220
Ambient conditions	24
Anesthetics	38
Application software installed on the computer	58
Application software user interface	12
Application software, installed on the computer	58, 62
Application software, main window	100
Application software, starting of	96
Application software, user interface	12
Applicator sterilizations, documentation	196
Applicators	127

Applicators, limited frequency of sterilization	127
Applicators, ordering as a set or single items	243
Applicators, selection of	163
Archive, not permanent	22
Archiving, non-permanently	32

B

Beam current, PRS control unit	248
Beam deflection adjustment	142
Beam deflection sensitivity	214
Beam deflection window	184
Beam deflection, graphic display	189
Before every use	29
Beryllium (Be)	64
Breast cancer	90

C

Cable connections	92
Cables - color coding	74, 93
Calibrating the IC ionization chamber	209
Calibrating the temperature sensor	223
Calibration intervals of X-ray radiation source	118
Calibration of XRS X-ray sources	209
Calibration reference	30
Calibration regulations	209
Calibration, screen mask	210
Calibrations	35, 210
Calibrations at the manufacturing plant	210
Cart	80
CD, burning to	102
Cellphones	252
Changing the line fuse	233
Checklist	178
Checklist for the OR	265
Class A	39
Class B (practice area)	25
Clean instruments immediately	238
Cleaning	39, 234, 235, 236
Cleaning and disinfecting the system cart	238
Cleaning and disinfection agents, unsuitable	238
Cleaning the touchscreen	237
Clinical consideration	36
Clinical experience	23
Clinical trials	23
Color coding	93

Comments, adding of	186
Component selection for treatment	164
Components to be checked	14
Components, cleaning and disinfecting	235
Computer level of the user terminal	106
Computer, in the monitor housing	58, 62
Conducting treatment	180
Configurations	24, 252
Connecting	92
Connections, various	58, 62
Continuous operation (standby)	53, 248
Control buttons for the main window	100
Control console	4
Control console display	54
Control lamps	34, 53
Control unit PRS 500	4
Controls, displays, connections	54
Conversion of the system	29
Core System, setting up	92
Count rate of ERM	184
Count rate of IRM	184
Current consumption	247

D

Damage and changes, obvious	86
Data export icons	106
Data import and export	106
Data import icons	106
Data loss	32
Data protection	31
Data protection responsibilities	31
Data, export	108, 110
Date and time synchronization	212
Date of manufacture	46
Deflecting coils, sensitivity of	214
Depth dose curve, measurement	76
Depth-dose rate	209
Designation in these Instructions for Use	4
Designations and abbreviations	4
Dilator Adjustment	226
Dilator adjustment using PAICH	226
Dimensions	249
Dimensions, PRS control unit	248
Dimensions, UNIDOS E	247
Disinfectants	238
Disinfectants, standard	238

Disinfecting surfaces	238
Disinfection	39, 234, 235, 236
Disinfection agents	238
Display for the setting parameters	53
Disposal	240
Disposal of components of the INTRABEAM Core System	41
Disposal of consumables	41
Dosage characteristics	246
Dosimeter UNIDOS E	72, 75
Drape	4
Drives for floppy disks and CD-ROM	58, 62
Dummy plug	56
Dynamic Offsets	142

E

Electrical features of the INTRABEAM Core System	37
Electrical requirements	86
Electrical safety	27, 30
Electromagnetic compatibility	25, 252
Electromagnetic interference	25, 39
EMC (electromagnetic compatibility)	39, 252
EMC requirements of IEC 60601-1-2	25, 39
Emergency backup power supply	27
EMERGENCY STOP	34
Emergency Stop button	188
Emergency stop button	54
Emission of an exact dose	73
End of the treatment session	190
Ending system operation	29
Enter your password	97, 173
Equipment Inventory	112
Equipment preparation	239
Equipment, capturing of	112
ERM check	161
Exporting data from the user terminal	106
External interlock	250
External radiation monitor ERM, checking	161

F

Field of application	22
Fire	38
Fuse change	233

G

Group, "XRS X-Ray Sources"	114
----------------------------------	-----

H

Handling	239
Handling and preparing the medical device	22
Hazard symbols	21
Heat development, excessive	181
High voltage generator	64

I

Icons for data export	106
Icons for data import	106
Importing data to the user terminal	106
Information security	31
Information symbols	21
Information, required additional	26
Initial calibrations	210
Input voltage	249
Input voltage, PRS control unit	248
Input voltage, UNIDOS E	247
Input window for comments	186
Instructions for Use, existing	26
Intended use	22
Interference with the system, unauthorized	25
Interlock connection, external	250
Interlock, external	34
Interrupting treatment	188
INTRABEAM cart	4, 26, 80, 264
INTRABEAM Core System	4, 12
INTRABEAM Core System, components for protection of	51
INTRABEAM Core System, disposal of components	41
INTRABEAM Core System, electrical features	37
INTRABEAM Core System, main components	50
INTRABEAM Core System, safety devices	33
INTRABEAM Core System, use of	30
INTRABEAM drape	4, 244
INTRABEAM Flat Applicator Set	27, 264
INTRABEAM needle applicator	26, 51, 244, 264
INTRABEAM radiation shields	4
INTRABEAM spherical applicators	26, 264
INTRABEAM stand	4
INTRABEAM Surface Applicator Set	27, 264
INTRABEAM System PRS 500 with XRS 4	4

INTRABEAM Water Phantom	27, 264
Intra-operative radiation therapy	22
Ionization chamber holder PAICH	70
IRM internal radiation monitor	65, 154
IRM linearity	218
IRM X-ray radiation monitor, internal	64
Isotropy	146, 153
Isotropy, adjust the	72
Isotropy, adjustment of	153
Isotropy, measurement	76

K

Kabel - Farbcodierung	74
Key to symbols	21

L

Label on the XRS 4 X-ray source	47
Layout diagram	89
Layout of equipment and personnel	89
Legal regulations	24
Liability and warranty	25
Line voltage required	86
Linearity of IRM internal radiation monitor, determination	218
Literature (supplementary documents)	266
Loss of warranty	25

M

Major performance characteristics	23
Malfunction	30
Managing XRS X-ray radiation sources	114
Manufacturer	46
Manufacturer's address	2
Measure the isotropy	72
Membrane keyboard with touchpad	58, 62
Memory space, insufficient	32
Menu bar	101
Menu, Tools	103
Minimal-invasive, highly focused treatment	22
Modifications, not permitted	25

N

NC32 INTRABEAM	4
NC32 INTRABEAM Floor Stand	27, 81, 264
Needle Applicator	40, 234

Needle applicator	244
Needle applicator, disposal of consumables	41
Non-sterile area	89
Normal use	22

O

Offline mode	98
Online mode	99
Operating mode	62
Operation	177
Optical Interlock system	30
OR equipment, checklist for	265
Order data for country-specific cables	241
Order number	244
Ordering data	241
Ordering data for INTRABEAM flat applicator	244
Ordering data for INTRABEAM surface applicator	244
Output Check (dose rate)	153
Oxygen	38

P

PAICH	70, 154
PAICH Output Check (dose rate)	153
Password rules	97
Patient ID, entering	163
Patient Information	168
Pause button	184
Pause mode	188
PDA photodiode array	72
PDA Source Check	146, 153
PDF files, creating and saving	102
Photo Diode Array (PDA)	72
Physician, radiotherapist	173
Physicist	173
Potential equalization	34, 56
Power connection	249
Power connection, PRS control unit	248
Power consumption	249
Power consumption, maximum	86
Power consumption, PRS control unit	248
Power switch	56, 58, 62
Powering up the INTRABEAM Core System	96
Preparations for use	83
Preparing and handling the medical device	22
Primary monitor	65

Principal components	53
Probe adjuster/ionization chamber holder PAICH	70
Probe adjustment using PAICH	138
Probe protector	68
Probe tip	64
Product Information	26
Product name	4
Protection components	66
Protective clothing	30
Protective ground contact	28
Protective ground contact, properly connected	86
PRS 500 control console	53
PRS 500 control console, rear panel of the	56
PRS 500 control console, technical data	247, 248
PRS 500 Control System	96
PRS 500 control unit	12, 53
PRS 500 Terminal	4
PRS 500 user terminal	27, 264
PRS500 control console	154
PTW company	2

Q

Qualification of users, advanced	22
Quality assurance devices	123
Quill icon	186

R

Radiation	36
Radiation characteristic, aspherical	65
Radiation distribution	246
Radiation dose, administered	190
Radiation dose, prescribed	190
Radiation exposure	247
Radiation monitor ERM, checking	161
Radiation monitor IRM, linearity determination	218
Radiation output, monitored	36
Radiation output, stability of the	246
Radiation shield	244
Radiation shielding	30
Radiation therapy treatment, system for	22
Radiation therapy, minimal-invasive, highly focused	22
Radio-frequency generating devices	25, 39
Radiotherapist, treating physician	173
Radiotherapy	85
Rated frequency	249

Rated frequency, PRS control unit	248
Rated frequency, UNIDOS E	247
Rating label	43
Ready button	184
Recalibration of XRS X-ray sources	120
Records, treatment and verification	53
Registered trademark of PTW	2
Regulatory information	260
Reports menu	133, 152, 160
Requirements for system operation	85
Responsibilities for data protection	31
Responsibilities of the owner/operator	24
Resterilization	234
Re-use	234
Review Treatment	199

S

Safety devices of the INTRABEAM Core System	33
Safety during tests and verification	35
Safety during treatment	29, 36
Safety measures	19
Saving the calibration file	120
Screen brightness, increase	62
Screen brightness, reduce	62
Screen mask, "Deflection Calibration"	215
Screen mask, "Enter Prescription"	163
Screen mask, "Equipment Inventory"	112
Screen mask, "Treatment"	181
Screen mask, "Verify Treatment Parameters"	174
Screen mask, Verify Treatment Parameters	173
Screen mask, "Deflection Factor History"	216
Selecting the XRS X-ray source	163
Sensor window	72
Service Identification Program (SIP)	227
Service requirements	227
Setting up the INTRABEAM Core System	92
Shields	4, 244
Shipment	239
Single-use material, disposal of consumables	41
SIP label	46
SIP label with SIP number	227
SIP number	46
Software and electronics	53
Solid water, water-equivalent material	72
Space requirement	85
Spare parts	242

Specialist staff, qualified	22
Specialist staff, trained and briefed	22
SQA reports, printing	133, 152, 160
SQA, System Quality Assurance	130
Staff, trained	85
Start button	184
Starting treatment	184
Start-up time, UNIDOS E	247
Sterile area	89
Sterilization	39, 234
Sterilization container, ordering as a single item.	243
Sterilization for interstitial treatment	39
Sterilization for tumor bed treatment	39, 234
Sterilization tray	78
Stop button	188
Stopping treatment	189
Straightness measurement, XRS probe and dilators	70
Supplementary documents	263
Surface temperature at the applicators	36
System calibration	210
System cart	4
System configuration	24, 252
System description	50
System for radiation therapy treatment	22
System Quality Assurance (SQA)	130
System Quality Assurance, SQA	131

T

Tab bar	105
Technical data for PRS 500 control console	247
Technical data of the PRS 500 control console	248
Technical data of the user terminal	249
Technical data of XRS 4 x-ray source	246
Tests and verification, safety during	35
Thermal energy, production of	36
Thermometer Calibration	223
Time and date synchronization	212
Tips for the owner/operator	24
Tools menu	103, 104
Touch function deactivation	62
Touchpad	58, 62
Touchscreen monitor	58, 62
Trademark	2
Training courses	85
Transport	28, 239
Transport and sterilization tray	78

Treatment facility	89
Treatment facility, layout	89
Treatment method, selection of	163
Treatment parameters, checking	173
Treatment parameters, display	199
Treatment parameters, entry of	163
Treatment record, personally signed	173
Treatment room	89
Treatment, interruption of	188
Treatment, safety during	29, 36
Treatment, stopping of	188, 189
Troubleshooting table	227
Troubleshooting	191
Tumor bed treatment, layout during	90
Tumor, direct radiation	23
Tumor, open surgical resection	23

U

UNIDOS E dosimeter	75, 154
UNIDOS E dosimeter from PTW	27, 264
UNIDOS E from PTW	4
UNIDOS E, start-up time	247
Use of the INTRABEAM Core System	30
User groups	97
User manuals, other applicable	264
User terminal	4, 12, 58
User terminal, technical data	249

V

Validity of verification	130
V-block guide	67
Verification components	50
Verification components, managing all	123
Verification devices	67, 70
Verify Treatment Parameters screen mask	173
Video instruction, video window	136
Visual system check	86

W

Warning and Information Labels	43
WARNING sign	45
Warranty and Liability	25
Water phantom	76
Water-equivalent material, solid water	72

Weight	249
Weight, PRS control unit	248
Weight, UNIDOS E	247
Workflow	94, 95
Working mode selection	98

X

X-block	154
X-block guide	66
X-ray exposure time	30
X-ray field	209
X-ray interlock, external	250
X-ray monitor, internal	36
X-ray probe	64
X-ray radiation	22
X-ray radiation display, external	251
X-ray radiation, defined, isotropic, constant	23
X-ray radiation, external display	250
X-ray source, status of	114, 127
X-ray source, technical data	246
XRS probe adjustment using PAICH	138
XRS X-ray radiation source	64
XRS X-ray radiation source, annual calibration	120
XRS X-ray source	4, 12, 154
XRS X-ray source, checking	76
XRS X-ray source, storage	66
XRS X-ray source, symbols on the	47
XRS X-ray sources	114

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Developer	
Name of the legal entity with transliteration	Carl Zeiss Meditec AG
Developer location address	Rudolf-Eber-Strasse 11, 73447 Oberkochen, Germany
Manufacturer	
Name of the legal entity with transliteration	Carl Zeiss Meditec AG
Legal entity location address	Goeschwitzer Strasse 51-52, 07745 Jena, Germany
Telephone numbers	Tel.: +49 36 41/220 439, Fax. +49 36 41/220 756
Legal entity e-mail	info@meditec.zeiss.com
Medical device manufacture place	Carl Zeiss Meditec AG, Rudolf-Eber-Strasse 11, 73447 Oberkochen, Germany

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Annex 2. Definitions

Application for registration of medical device	Operational documentation, extract from the manufacturer's technical documentation
I. INTRABEAM Intraoperative Radiotherapy System, version INTRABEAM PRS 500 Consisting of:	INTRABEAM system, INTRABEAM PRS 500, system
1. X-ray source (XRS 4).	XRS source, X-ray source XRS, radiation source XRS 4, XRS 4
2. V-Block.	"V-block" guide
3. X-Block.	X-block
4. Control console PRS 500.	Control console PRS 500
5. Terminal PRS 500 (Computer with touchscreen and INTRABEAM Software).	User terminal PRS 500, terminal PRS 500
6. Probe Adjuster Ion Chamber Holder (PAICH).	Probe centering device / holder PAICH, PAICH (probe centering device and ionization chamber holder)
7. Ionization chamber / Ionization Chamber Shield.	Ionization chamber IC, Ionization chamber "IC shield"
8. Photo diode array (PDA).	Photodiode array PDA.
9. Electrometer (UNIDOS E) manufactured by PTW Freiburg Physikalisch-Technische Werkstätten Dr. Pöchlau GmbH, production site located at: Lörracher Str. 7, 79115 FREIBURG, GERMANY.	Dosimeter UNIDOS E, Electrometer UNIDOS E, UNIDOS E
10. INTRABEAM Water Phantom.	Water Phantom, INTRABEAM Water Phantom
11. Isolating Transformer.	---
12. INTRABEAM Cart.	INTRABEAM Workplace, Cart INTRABEAM
13. NC32 INTRABEAM Floor Stand.	Floor Stand NC32 INTRABEAM, Floor Stand
14. Applicators for INTRABEAM Intraoperative Radiotherapy System in sets and on a standalone basis, with accessories P3H 2014/2108 of 01.11. 2018 r. - INTRABEAM Spherical Applicator 1.5 cm (set of INTRABEAM Spherical Applicator 1.5 cm, calibration disc) - INTRABEAM Spherical Applicator 2.0 cm (set of INTRABEAM Spherical Applicator 2.0 cm, calibration disc)	INTRABEAM Spherical Applicator, Spherical Applicator

- INTRABEAM Spherical Applicator 2.5 cm (set of INTRABEAM Spherical Applicator 2.5 cm, calibration disc)
- INTRABEAM Spherical Applicator 3.0 cm (set of INTRABEAM Spherical Applicator 3.0 cm, calibration disc)
- INTRABEAM Spherical Applicator 3.5 cm (set of INTRABEAM Spherical Applicator 3.5 cm, calibration disc)
- INTRABEAM Spherical Applicator 4.0 cm (set of INTRABEAM Spherical Applicator 4.0 cm, calibration disc)
- INTRABEAM Spherical Applicator 4.5 cm (set of INTRABEAM Spherical Applicator 4.5 cm, calibration disc)
- INTRABEAM Spherical Applicator 5.0 cm (set of INTRABEAM Spherical Applicator 5.0 cm, calibration disc)
- INTRABEAM Flat Applicator 1.0 cm (set of INTRABEAM Flat Applicator 1.0 cm, INTRABEAM Position Marker 1.0cm, INTRABEAM Lumen Plug, calibration disc)
- INTRABEAM Flat Applicator 2.0 cm (set of INTRABEAM Flat Applicator 2.0 cm, INTRABEAM Position Marker 2.0cm, INTRABEAM Lumen Plug, calibration disc)
- INTRABEAM Flat Applicator 3.0 cm (set of INTRABEAM Flat Applicator 3.0 cm, INTRABEAM Position Marker 3.0cm, INTRABEAM Lumen Plug, calibration disc)
- INTRABEAM Flat Applicator 4.0 cm (set of INTRABEAM Flat Applicator 4.0 cm, INTRABEAM Position Marker 4.0cm, INTRABEAM Lumen Plug, calibration disc)
- INTRABEAM Flat Applicator 5.0 cm (set of INTRABEAM Flat Applicator 5.0 cm, INTRABEAM Position Marker 5.0cm, INTRABEAM Lumen Plug, calibration disc)
- INTRABEAM Flat Applicator 6.0 cm (set of INTRABEAM Flat Applicator 6.0 cm, INTRABEAM Position Marker 6.0cm, INTRABEAM Lumen Plug, calibration disc)
- INTRABEAM Surface Applicator 1.0 cm (set of INTRABEAM Surface Applicator 1.0 cm, INTRABEAM Position Marker 1.0cm, INTRABEAM Lumen Plug, calibration disc)
- INTRABEAM Surface Applicator 2.0 cm (set of INTRABEAM Surface Applicator 2.0 cm, INTRABEAM Position Marker 2.0cm, INTRABEAM Lumen Plug, calibration disc)
- INTRABEAM Surface Applicator 3.0 cm (set of INTRABEAM Surface Applicator 3.0 cm,

INTRABEAM Flat Applicator, Flat Applicator

INTRABEAM Surface Applicator, Surface Applicator

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9. Photo diode array (PDA)	Photodiode array PDA.
10. Electrometer (UNIDOS E) manufactured by PTW Freiburg Physikalisch-Technische Werkstätten Dr. Pöchlauer GmbH, production site located at: Lörracher Str. 7, 79115 FREIBURG, GERMANY.	Dosimeter UNIDOS E, Electrometer UNIDOS E, UNIDOS E
11. INTRABEAM Water Phantom	Water Phantom, INTRABEAM Water Phantom
12. Isolating Transformer	---
13. INTRABEAM Workplace	INTRABEAM Workplace, INTRABEAM Cart
14. NC32 INTRABEAM Floor Stand	Floor Stand NC32 INTRABEAM, Floor Stand
15. Applicators for INTRABEAM Intraoperative Radiotherapy System in sets and on a standalone basis, with accessories P3H 2014/2108 of 01.11. 2018 r. - INTRABEAM Spherical Applicator 1.5 cm (set of INTRABEAM Spherical Applicator 1.5 cm, calibration disc) - INTRABEAM Spherical Applicator 2.0 cm (set of INTRABEAM Spherical Applicator 2.0 cm, calibration disc) - INTRABEAM Spherical Applicator 2.5 cm (set of INTRABEAM Spherical Applicator 2.5 cm, calibration disc) - INTRABEAM Spherical Applicator 3.0 cm (set of INTRABEAM Spherical Applicator 3.0 cm, calibration disc) - INTRABEAM Spherical Applicator 3.5 cm (set of INTRABEAM Spherical Applicator 3.5 cm, calibration disc) - INTRABEAM Spherical Applicator 4.0 cm (set of INTRABEAM Spherical Applicator 4.0 cm, calibration disc) - INTRABEAM Spherical Applicator 4.5 cm (set of INTRABEAM Spherical Applicator 4.5 cm, calibration disc) - INTRABEAM Spherical Applicator 5.0 cm (set of INTRABEAM Spherical Applicator 5.0 cm, calibration disc) - INTRABEAM Flat Applicator 1.0 cm (set of INTRABEAM Flat Applicator 1.0 cm, INTRABEAM Position Marker 1.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Flat Applicator 2.0 cm (set of INTRABEAM Flat Applicator 2.0 cm, INTRABEAM Position Marker 2.0cm, INTRABEAM Lumen Plug, calibration disc)	INTRABEAM Spherical Applicator, Spherical Applicator INTRABEAM Flat Applicator, Flat Applicator

- INTRABEAM Flat Applicator 3.0 cm (set of INTRABEAM Flat Applicator 3.0 cm, INTRABEAM Position Marker 3.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Flat Applicator 4.0 cm (set of INTRABEAM Flat Applicator 4.0 cm, INTRABEAM Position Marker 4.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Flat Applicator 5.0 cm (set of INTRABEAM Flat Applicator 5.0 cm, INTRABEAM Position Marker 5.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Flat Applicator 6.0 cm (set of INTRABEAM Flat Applicator 6.0 cm, INTRABEAM Position Marker 6.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Needle Applicator (set of INTRABEAM Needle Applicator, guide shaft tapered (3 pcs), guide shaft beveled (3pcs)) - INTRABEAM Surface Applicator 1.0 cm (set of INTRABEAM Surface Applicator 1.0 cm, INTRABEAM Position Marker 1.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Surface Applicator 2.0 cm (set of INTRABEAM Surface Applicator 2.0 cm, INTRABEAM Position Marker 2.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Surface Applicator 3.0 cm (set of INTRABEAM Surface Applicator 3.0 cm, INTRABEAM Position Marker 3.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Surface Applicator 4.0 cm (set of INTRABEAM Surface Applicator 4.0 cm, INTRABEAM Position Marker 4.0cm, INTRABEAM Lumen Plug, calibration disc)	INTRABEAM Needle Applicator, Needle Applicator INTRABEAM Surface Applicator, Surface Applicator
16. Container for sterilizing INTRABEAM Spherical Applicators (including re-usable filter for container)	Sterile container
17. Transport case	-
18. Dustcover	-
19. INTRABEAM DRAPE, place of production Microtek Dominicana, S.A., production site located at: Zona Franca II, La Romana 22000, Dominican Republic	-
20. Portable data storage device	
21. XRS Cable (cable for connection of INTRABEAM 600 and XRS radiation source)	cable XRS
22. Network/mains cable	
23. PAICH/PDA Cable	-
24. Radiation therapy treatment planning simulation system: Radiance software, manufactured by GMV Soluciones Globales	program «radiance», radiance

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
Reiner Hammel

Certified Copy


I hereby certify that the above copy is a true copy of the original, which has been produced by me.

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2019-11-14

ZEISS

- About these Instructions for Use* The Instructions for Use are part of the delivery package.
- Please read them carefully before using the device.
 - Keep them at the place where the system is used.
 - Keep them for the entire service life of the system.
 - Pass them on to any subsequent owner or user of the system.
- Scope* The present Instructions for Use apply to INTRABEAM 600 with REF number 6405.
- Trademark* INTRABEAM is a registered trademark of Carl Zeiss Meditec AG or of another company of the ZEISS Group in Germany and/or other countries.
- Trademark* UNIDOS®
is a registered trademark of the PTW company.
- radiance®
is a registered trademark of the GMV company, Spain.
- Manufacturer* Carl Zeiss Meditec AG
Goeschwitzer Strasse 51-52
07745 Jena
Germany
- Email: surgical@meditec.zeiss.com
Internet: www.zeiss.com/med
- Original language* These Instructions for Use were originally drafted and published in German.
- Copyright* © Carl Zeiss Meditec AG 2019
All rights reserved. Reprints and reproductions - including excerpts - are not permitted. However, this user manual may be digitized and archived for internal use. Access to this data by third parties is not permitted and must be prevented.
- Change note* Subject to change in design and scope of delivery and as a result of ongoing technical development.

Printed in Germany.

Designations and abbreviations

Several official (English) product names have been abbreviated in the respective national language in order to enhance the clarity and readability of these Instructions for Use. In cases where variants of individual devices and products are explained, these are clearly identified in the text.

Name of product (English)	Designation in these Instructions for Use	Explanation
INTRABEAM 600	INTRABEAM 600	INTRABEAM 600 refers to the complete product, workplace and software.
INTRABEAM Workplace	Workplace	The workplace is the operating and safe-keeping station of the INTRABEAM 600. It is equipped with the CC 600, a computer, a monitor, a dosimeter and the stored components XRS 4, PAICH, PDA.
XRS 4 X-ray source	XRS 4	XRS (X-Ray Source), functional principle: linear accelerator; The XRS 4 is stored in the workplace.
Control Console 600	CC 600	Control unit and power supply for the XRS 4
Computer	Computer	Permanently installed in the workplace, access to CD/DVD drive is possible
Touchscreen monitor	Monitor	Current size: 27" Permanently installed in the workplace
UNIDOS E dosimeter	Dosimeter	Currently: UNIDOS E dosimeter manufactured by PTW, permanently installed in the workplace
PAICH	PAICH	Probe adjuster and ionization chamber holder, stored in the workplace
PDA	PDA	Photo Diode Array , stored in the workplace
IC shield	Ionization chamber	Ionization chamber with shield, stored in the workplace
NC32 INTRABEAM Floor Stand	Stand	Accessory for the INTRABEAM 600, currently: NC32 INTRABEAM Floor Stand
INTRABEAM Drape	Drape	Sterile covering for the INTRABEAM Stand

Name of product (English)	Designation in these instructions for Use	Explanation
INTRABEAM Spherical Applicator	Spherical applicator	Spherical applicator
INTRABEAM Surface Applicator	Surface applicator	Surface applicator
INTRABEAM Flat Applicator	Flat applicator	Flat applicator
INTRABEAM Needle Applicator	Needle applicator	Needle-shaped applicator, sterile disposable product

Abbreviations in these Instructions for Use

EU	European Union
IC	Ionization chamber
IRM	Interner Röntgenstrahlen-Monitor (internal X-ray monitor)
PAICH	Probe adjuster and ionization chamber holder
PDA	Photodiode array
QA	Quality assurance
EOT	End of treatment
RTD	Resistance temperature device
EMC	Electromagnetic compatibility
PTB	Physikalisch Technische Bundesanstalt (Federal Institute of Physics and Technology)
NIST	U.S. National Institute of Standards and Technology
USB	Universal Serial Bus (serial bus system of the computer)
ABC error	Automatic beam centering error
HVL	Half-value layer

Contents

Designations and abbreviations	3
Abbreviations in these Instructions for Use	4

Functions at a glance

INTRABEAM 600	14
Components for inspection (in the workplace)	16
Accessories for inspection(not in the workplace)	16
Accessories for the INTRABEAM 600	18

Safety measures

Key to symbols	23
Hazard symbols	23
Information symbols	23

Target group	24
--------------------	----

Field of application	24
Intended use	24
Indications for use	24
Major performance characteristics	25
Normal use	25

Information for operators and users	27
---	----

Use of the INTRABEAM 600	34
--------------------------------	----

Data protection and information security	35
Archiving of data	36
Data carriers must be free from viruses.	36

Safety devices of the INTRABEAM 600	37
Optical interlock system	37
EMERGENCY STOP, control lamps and external interlock	38
Calibrations	39
Safety during tests and verification	39
Safety during treatment	39
Electrics of the INTRABEAM 600	40
Precautions regarding EMC	42
Fire, anesthetics and oxygen	43
 Cleaning, disinfecting, sterilizing	 44
 INTRABEAM 600 product identification	 46
Warning and information labels on the INTRABEAM 600 Workplace	46
Inscriptions, labels and symbols on the components	56
Label on battery compartment of UNIDOS E	62

System overview

System description	64
Overview	64
 INTRABEAM Workplace	 66
Connectors on the workplace	68
 Principal components	 69
Control Console 600 (CC 600)	69
Computer of the INTRABEAM 600	72
XRS 4 with IRM internal X-ray monitor	74
 Protection components	 77
Probe protector	77
V-guide	78

Verification devices	80
PAICH probe adjuster and ionization chamber holder	80
Ionization chamber "IC shield"	81
PDA (photo diode array)	83
UNIDOS E dosimeter	84
 Accessories for your INTRABEAM 600	 86
NC32 INTRABEAM Floor Stand	86
Applicators	88
INTRABEAM Water Phantom	89

Preparations for use

Requirements for system operation	93
 Preparing the workplace for use	 94
Transporting the workplace	94
Connecting the workplace to the power supply	97
Visual system check	97
 Spatial arrangement of the treatment facility	 99
 Connecting	 102
 Workflow	 104
 Structure of the software user interface	 106
Tab bar	107
Menu bar	108
Buttons	108

Startup of the INTRABEAM 600	111
Switching the workplace on	111
Starting the computer	112
INTRABEAM Start Manager (start window)	113
1 <INTRABEAM> button	114
Select working status	116
2 "radiance" software	117
3 "LINAC Configuration tool" software	118
4 Importing calibration data	119
5 Data Management	120
 Menus in the menu bar	 124
1 <File> menu	124
2 <Tools> menu	125
3 <QA Reports> menu	133
4 Menu item <Patient Database>	133
5 Menu item <About>	138
 Equipment inventory	 140
Equipment inventory, "XRS linear accelerator" group	142
Recalibrating the XRS 4	146
Equipment inventory, "Equipment" group	148
Equipment inventory, "Applicators" group	152
 Quality assurance	 158
Performing quality assurance	160
Video instruction, video window	166
XRS tip adjustment	168
Dynamic Offsets (adjusting the beam deflection)	173
Isotropy check	176
Output check of dose rate	182

Treatment prescription	190
Selecting a patient	190
Marking the treatment site	192
Entering treatment parameters	194
Planning fractional treatment	202
Treatment approval	205

Operation

Checklist for the INTRABEAM 600	210
Conducting treatment	212
"Treatment monitoring" window	216
Starting the treatment	218
Interrupting the treatment (Pause mode)	222
Canceling a treatment (Stop)	224
Treatment terminated	226
Troubleshooting	227
Review treatment, displaying treatment parameters	231
Displaying the treatment record	242
Generating a treatment record	244

Care and maintenance

Calibration Instructions	254
Calibration	257
Determining the beam centering	258
Determining the linearity of the IRM internal X-ray monitor	262
Calibrating the air pressure sensor	264
Calibrating the temperature sensor	267

Safety check	270
Workplace safety check	270
Troubleshooting table	271
Changing the line fuse of the workplace	277
Handling, transport and shipment	278
Handling components of the INTRABEAM 600	278
Care of the INTRABEAM 600	279
Handling components	281
Cleaning components	282
Disinfecting components	284
Maintenance intervals	285
Disposal	288

System data

Order data	290
Accessories	291
Spare parts	292
Technical data	293
Technical data of the INTRABEAM Workplace	293
Technical data of isolating transformer	294
Technical data of the XRS 4	295
Technical data for UNIDOS E	296
Technical data for Control Console 600 (CC 600)	297
Technical data of computer	297
Ethernet connection (optional),	298
Data carriers must be virus-free.	299

External interlock and X-radiation indicator	300
External X-ray display	302
EMC (electromagnetic compatibility)	303
EMC - Electromagnetic compatibility as per IEC 60601-1-2: 2004/2007 (2nd/3rd Edition)	305
EMC - Electromagnetic compatibility as per IEC 60601-1-2: 2014 (4th Edition)	309
Environmental conditions	311
Regulatory information	312

Supplementary documents

Applicable user manuals	314
Literature on the INTRABEAM 600	315
Checklist for OR equipment	316

Index


This image shows a single sheet of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page. There is no handwriting or other markings on the paper.

Functions at a glance

INTRABEAM 600	14
Components for inspection (in the workplace).....	16
Accessories for inspection (not in the workplace).....	16
Accessories for the INTRABEAM 600.....	18

INTRABEAM 600

INTRABEAM Workplace		Page 69
<u>1</u>	Control Console 600 (CC 600)	Page 70
<u>2</u>	Connection for ionization chamber	Page 70
<u>3</u>	Drawer i.a. with Instructions for Use	Page 70
<u>4</u>	Shielding with ionization chamber	Page 81
<u>5</u>	PAICH	Page 80
<u>6</u>	PDA	Page 83
<u>7</u>	XRS 4 X-ray source with integrated IRM X-ray monitor and probe protector (7.1)	Page 102
<u>8</u>	Cable for PDA/PAICH	Page 102
<u>9</u>	Cable for XRS	Page 102
<u>10</u>	EMERGENCY STOP switch on the control unit	Page 69
<u>11</u>	Drawer with keyboard + mouse	Page 72
<u>12</u>	Computer	Page 102
<u>13</u>	UNIDOS E Dosimeter	Page 102
Software user interface		Page 106
Tab for calling up important function groups:		
<u>14</u>	<Equipment Inventory> tab	Page 140
<u>15</u>	<Quality Assurance> tab	Page 158
<u>16</u>	<Treatment Prescription> tab	Page 190
<u>17</u>	<Verify Treatment Parameters> tab	Page 205
<u>18</u>	<Treatment> tab	Page 231
<u>19</u>	<Review Treatment> tab	Page 231
<u>20</u>	<Calibration> tab	Page 257
Controls		
<u>21</u>	Touchscreen monitor	
<u>22</u>	Mouse	
<u>23</u>	PC keyboard (membrane keypad)	
<u>24</u>	USB port	


Components for inspection (in the workplace)

1	V-guide The V-guide is embedded in the working surface of the workplace and is removable.	Page 78
2	PDA photo diode array, stored in the workplace	Page 83
3	PAICH probe adjuster and ionization holder, stored in the workplace	Page 80
4	Ionization chamber with shield, stored in the workplace	Page 81
5	UNIDOS E dosimeter, integrated in the workplace	Page 84

Accessories for inspection (not in the workplace)

6	INTRABEAM Water Phantom, option	Page 89
---	------------------------------------	---------


Accessories for the INTRABEAM 600


As the products are not licensed in all countries, please contact the local distributors with regard to the licenses.

Information on the accessories of the INTRABEAM system can be found in the corresponding separate Instructions for Use. Overview, see Page 30.

- | | |
|---|--|
| 1 | INTRABEAM Spherical Applicator
(The illustration shows the applicator set in the sieve container) |
| 2 | INTRABEAM Flat Applicator |
| 3 | INTRABEAM Surface Applicator |
| 4 | INTRABEAM Needle Applicator |
| 5 | Radiation shields |
| 6 | NC32 INTRABEAM Floor Stand |
| 7 | INTRABEAM Drapes (sterile covers) |

Blank lined paper.

Safety measures

Key to symbols	23
Hazard symbols	23
Information symbols	23
Target group	24
Field of application	24
Intended use	24
Indications for use	24
Major performance characteristics	25
Normal use	25
Information for operators and users	27
Use of the INTRABEAM 600	34
Data protection and information security	35
Archiving of data	36
Data carriers must be free from viruses.	36
Safety devices of the INTRABEAM 600	37
Optical interlock system	37
EMERGENCY STOP, control lamps and external interlock	38
Calibrations	39
Safety during tests and verification	39
Safety during treatment	39
Electrics of the INTRABEAM 600	40
Precautions regarding EMC	42
Fire, anesthetics and oxygen	43

Cleaning, disinfecting, sterilizing	44
INTRABEAM 600 product identification	46
Warning and information labels on the INTRABEAM 600 Workplace	46
Inscriptions, labels and symbols on the components	56
Label on battery compartment of UNIDOS E	62

Key to symbols

We would like to inform you about safety aspects which must be observed when handling this device. This chapter contains a summary of the most important information concerning matters relevant to instrument safety.

Hazard symbols

The following safety information has been incorporated into the Instructions for Use. Please note this information and be particularly careful in these cases.

WARNING

Indicates a hazard which **can lead to death or severe injury** if it is not prevented.

CAUTION

Indicates a hazard which can cause **minor to moderate injury** if it is not prevented.

NOTE

Indicates a hazard which can cause **damage to material** if it is not prevented.

Information symbols

The following information symbols are used in these Instructions for Use:

- Listing
- Prompt for action
- ✓ Prerequisite for an action
- Result of an action

Additional information and tips

Target group

These Instructions for Use are directed towards qualified specialists who have been trained and briefed in preparing and handling this medical device using the Instructions for Use.

Certain fields of application of the INTRABEAM 600 may require the user to obtain additional qualifications in the corresponding areas of expertise.

Field of application

Intended use

The INTRABEAM 600 is intended for use in radiotherapy treatments.

Indications for use

The INTRABEAM 600 is indicated for radiation therapy treatments.

- The INTRABEAM Spherical Applicators are indicated for use with the INTRABEAM 600 to deliver a prescribed dose of radiation to the treatment margin or tumor bed during intracavity and intraoperative radiotherapy treatments.
- The INTRABEAM Needle Applicator (comprising of Needle Applicator and guide shafts) is intended for use in combination with the INTRABEAM 600 to intraoperatively administer radiation to tissue including irradiation of intracranial tumors.

- The INTRABEAM Flat Applicator Sets (each comprising of INTRABEAM Flat Applicator, INTRABEAM Position Marker and INTRABEAM Lumen Plug) is intended to supply a specified radiation dose during applications in combination with the INTRABEAM 600,
 - during intraoperative radiotherapy, on a surgically exposed surface or in a tumor bed,
 - during treatment of tumors on the body surface.

The INTRABEAM Flat Applicator Sets is designed to deliver a flat radiation field at a distance of 5 mm (in water) from its circular application surface.

- The INTRABEAM Surface Applicator Sets (each comprising of INTRABEAM Surface Applicator, INTRABEAM Position Marker and INTRABEAM Lumen Plug) is intended to supply a specified radiation dose during applications in combination with the INTRABEAM 600,
 - during intraoperative radiotherapy, on a surgically exposed surface or in a tumor bed.
 - during treatment of tumors on the body surface.

The INTRABEAM Surface Applicator Sets is designed to deliver a flat radiation field directly at the applicator's surface.

Major performance characteristics

The major performance characteristic of the INTRABEAM 600 is the delivery of defined, isotropic and constant X-ray radiation.

Normal use

The INTRABEAM 600 can be used for minimally invasive, highly-focused treatment of selected lesions and for intraoperative radiotherapy of a tumor bed throughout the entire body.

The core piece of the INTRABEAM 600 is the XRS 4 miniature X-ray source. The XRS 4 is a linear accelerator used to generate X-rays.

The XRS 4 is equipped with an X-ray probe which is 10 cm long and has a diameter of 3.2 mm. This X-ray probe releases low-energy radiation in high doses. To do this, electrons are accelerated along the X-ray probe with a maximum acceleration voltage of 50 kV (linear accelerator). They strike a

target on the tip of the X-ray probe, thus generating X-ray radiation. The effectively emitted X-ray energy equals approx. 20keV or less. The dose decreases quickly due to the low energy of the X-rays in soft tissue (approx. $1/r^3$). The treated lesion can thus shield the healthy tissue against damage.

Various applicators are available for connection to the probe tip of the XRS 4 for placement on a lesion for local irradiation or insertion in the bed of a removed tumor.

The INTRABEAM 600 can be configured so as to enable the irradiation of defined lesions using interstitial, intraoperative and intracavitary procedures and the treatment of surfaces.

The NC32 INTRABEAM Floor Stand is available to you for reliable positioning of the XRS 4 with the applicator mounted on the surface. The XRS 4 and the stand arm can be sterilely covered using the INTRABEAM Drape.

Thanks to the movable INTRABEAM Workplace, the INTRABEAM 600 is a mobile workstation. The workplace is the control and safekeeping station of the INTRABEAM 600. The computer, CC 600 and dosimeter as well as all components for quality assurance and accessories are housed in the workplace.

The INTRABEAM Water Phantom, which you can use to check the depth dose curve and isotropy of your XRS 4 yourself, is available to you as an external quality assurance component.

ZEISS supplies routines for quality assurance and systematically recurring calibration.

The INTRABEAM 600 is not intended to be used for archiving data.

Literature on the INTRABEAM 600

Current clinical, physical and scientific publications on the INTRABEAM 600 and the respective fields of application can be found under www.zeiss.com/radiotherapy

Information for operators and users

Obligations of the operator

The INTRABEAM 600 may be used only by qualified specialists who have been trained in the operation of the INTRABEAM 600. No one else may operate or use the INTRABEAM 600.

Every facility that operates an INTRABEAM 600 must appoint suitable and qualified specialists who have been trained in operation of the INTRABEAM 600 and are responsible for the tasks associated with the INTRABEAM 600.

The manufacturer, Carl Zeiss Meditec AG, offers training courses on the INTRABEAM 600.

Only physicians, other authorized medical professionals or other persons designated by them may purchase the INTRABEAM 600.

- Keep the Instructions for Use readily available to the operating personnel at all times.

Legal regulations

- Observe the legal regulations for accident prevention and occupational health and safety applicable in the country concerned.

Medical electrical systems

Any additional equipment connected to medical electrical devices must demonstrably comply with the applicable standards and directives (e.g. IEC 60950-1:2005, modified + Cor.:2006 + A1:2009, modified + A1:2009/Cor.:2012 + A2:2013, modified for data processing equipment).

Furthermore, any configurations must comply with normative requirements for electronic medical devices (see IEC 60601-1:2005 + Cor. :2006 + Cor. :2007 + A1:2012).

Anyone connecting additional equipment to medical electrical devices is a system configurer and as such is responsible for compliance of the system with the normative requirements for medical electrical systems.

Please note that local legislation takes priority over the above-mentioned normative requirements. If you have any questions, please contact your local dealer or ZEISS Service.

- Only use original accessories or consumables approved by Carl Zeiss. If you wish to use other accessories ensure that Carl Zeiss or the manufacturer of the accessories have verified and confirmed that these accessories meet the respective safety standards and can be used without risk.

- Environmental conditions*
- Do not operate the INTRABEAM 600:
 - in potentially explosive atmospheres,
 - less than 25 cm away from any flammable anesthetics or volatile solvents such as alcohol, benzine or similar substances.
 - Do not use or store components of the INTRABEAM 600 in damp rooms.
 - Prevent dripping, surging or spraying water in the vicinity of components of the INTRABEAM 600.
 - Please ensure that no liquid is permitted to penetrate into the components of the INTRABEAM system.
 - Ensure that the installation conditions and the use of the device meet surgical requirements:
 - Low vibration
 - Clean environment
 - Avoid extreme mechanical stress
 - Observe the environmental conditions for operation, transportation and storage in the "Device details" section, Page 311.

Electromagnetic compatibility (EMC)

- Please note the information on electromagnetic compatibility (EMC) in the section on "Device data", Page 303.

The electrically operated components of the INTRABEAM 600 comply with

 - the EMC requirements of IEC 60601-1-2.
 - for class B requirements for electromagnetic interference (practices).
- It cannot be ruled out, however, that interferences may occur in RF receivers (e.g. TV or radio sets) in the vicinity. Should interferences be detected, please notify ZEISS Service.
- When operating the device please observe the following precautions regarding EMC:
 - Use only accessories which have been approved for this system by ZEISS.
 - Do not use any portable or mobile RF communication equipment in the vicinity of the INTRABEAM 600, as the possibility that this may impair the device's function cannot be excluded.

- Warranty and liability** The warranty and liability depend on the contractually specified conditions. The manufacturer is not liable for damage caused by unauthorized interference with components of the INTRABEAM 600. Furthermore, this will forfeit any rights to claim under warranty.
- Do not modify the components of the INTRABEAM 600 without the manufacturer's express permission.
 - If a component of the INTRABEAM 600 is modified after consultation with the manufacturer, suitable inspections and testing must be completed to ensure subsequent safe use.
 - Modifications and repairs of the INTRABEAM 600 or any equipment operated together with this device may only be performed by ZEISS Service or other specialists authorized by Carl Zeiss.
- Requirements for the user** Correct operation of the components of the INTRABEAM 600 is essential for safe operation of the INTRABEAM 600. Therefore, please thoroughly familiarize yourself with the content of these instructions for Use before starting up the INTRABEAM 600. Please also observe the Instructions for Use of any additional equipment and the accessory components. Supplementary information is available from our Service department or from authorized distributors.
- Use the INTRABEAM 600 only within the limits of its intended purpose.
- Secure network** The user (IT manager) is responsible for ensuring that the network is connected to the system via a suitable security design (e.g. firewall) to protect it against unauthorized access from outside.
- USB media, virus-free** The user is responsible for ensuring that the USB medium used for data exchange is virus-free.
- Note the symbols and labels that appear on the packaging and components of the INTRABEAM 600.

Required additional information

The present Instructions for Use contain detailed explanations on the use of the INTRABEAM 600 and information on the accessories.

For more information on the use of the INTRABEAM 600, please refer to the Instructions for Use or product information enclosed with the individual components and the accessories.

- Also thoroughly read through the following Instructions for Use before using the INTRABEAM 600:
 - NC32 INTRABEAM Floor Stand
ZEISS order number: 1096-945
Instructions for Use, ZEISS: G-30-1476
 - INTRABEAM Drape
ZEISS order number: 32090-0000-000
Instructions for Use, ZEISS: G-30-1526
 - Radiation shield, flat
ZEISS order number: 304534-6000-595
Instructions for Use, ZEISS: G-30-1890
 - INTRABEAM Spherical Applicator 1.5 cm, 2.0 cm, 2.5 cm, 3.0 cm, 3.5 cm, 4.0 cm, 4.5 cm, 5.0 cm and Set
ZEISS order numbers: 304534-6000-557, -558, -559, -560, -561, -562, -563, -564 and -348 (set)
Instructions for Use, ZEISS: G-30-1478
 - INTRABEAM Needle Applicator
ZEISS order number: 304534-3001-200
Instructions for Use, ZEISS: G-30-1723
 - INTRABEAM Surface Applicator Set 1.0 cm, 2.0 cm, 3.0 cm and 4.0 cm
ZEISS order numbers: 304534-3001-530, -531, -532, -533
Instructions for Use, ZEISS: G-30-1858
 - INTRABEAM Flat Applicator Set 1.0 cm, 2.0 cm, 3.0 cm, 4.0 cm, 5.0 cm and 6.0 cm
ZEISS order numbers: 304534-3001-565, -566, -567, -568, -569 and -570
Instructions for Use, ZEISS: G-30-1880
 - UNIDOS E dosimeter from PTW, Type 10010
PTW Instructions for Use, English: D545.131.00/06 en
PTW Instructions for Use, German: D545.151.00/06 de
 - INTRABEAM Water Phantom
ZEISS order number: 304534-7000-600
Instructions for Use, ZEISS: G-30-1767.

- Electrical safety**
- Always switch off the device before connecting it to or disconnecting it from the power grid, if it will not be used for a prolonged period of time, or for cleaning its surface.
 - Switch the device off at the power switch or disconnect the power plug if you notice any smoke, sparks or unusual noises. Do not use the device before it has been repaired by our Service.

A potential equalization connection can be found on the connection panel (see IEC 60601-1:2005 + Cor. :2006 + Cor. :2007 + A1:2012). This enables other active devices to be connected to the same ground potential, or a redundant connection to be made to the protective ground.

- Do not use any multiple sockets!
- Never open the device! The device contains freely accessible live components. If you remove the housing, you run the risk of electric shock.
- Do not touch the system if your body is electrostatically charged and the system is not grounded.
- Connect the device to the "emergency backup power supply" according to the country-specific regulations or guidelines.

- Transport & service**
- Over longer distances (e.g. removal, return for repair, etc.) the system must always be transported in the original packaging or in special return packaging. For details, please contact your dealer or the Carl Zeiss service team.
 - In order to prevent any impairment of the device's safety due to age, wear, etc., the user must ensure that the device is subjected to the necessary safety checks (see section "Care and Maintenance").

- ZEISS Service**
- Have the device inspected regularly by ZEISS Service. This device is a high quality technical product. We recommend this in order to ensure its trouble-free and safe functioning. Carl Zeiss Meditec AG offers corresponding service contracts.
 - If a failure occurs that you cannot correct according to the table "Troubleshooting", attach a sign to the device stating that it is out of order and contact ZEISS Service.

- Info to ZEISS Service* ZEISS Service have the service manual "SM-30-8023 INTRABEAM System" at their disposal for performing service. This service manual is designed only for internal use by ZEISS Service.
- Service Identification Program (SIP)* If your device or system has been entered in the ZEISS "Service Identification Program" (SIP), it will bear an SIP label with an SIP number, see Page 50. The SIP number is a unique identification number assigned to your device or system. ZEISS keeps a product file under this SIP number.
- For any service activities required, please turn to your local contact person of the ZEISS sales organization in your country. The telephone number, fax number and e-mail address can be found on the SIP label.
- Prerequisites for operation* Our service representative or an expert authorized by us will install the system. Please ensure that the following requirements for operation are maintained:
- All mechanical connections (details in the user's manual) which are relevant to safety are properly connected and screw connections tightened.
 - All cables and plugs are in good working condition.
 - The voltage setting on the instrument conforms to the rated voltage of the line supply on site.
 - The instrument is plugged into a power outlet which has a properly connected protective ground contact.
 - The power cord being used is the one designed for use with this instrument.
 - Do not force cable connections. If the male and female parts do not readily connect, make sure that they are appropriate for one another. If any of the connectors are damaged, have our service representative repair them.
 - Do not pull on the cables. The plug connections are self-locking. Slide back the sleeve to separate the connection.
- Before every use and after every system change*
- Ensure that all specified "requirements for operation" are met.
 - Work through the checklist, see Page 210.
 - Note all symbols, notes and labels on the system.

Safety during treatment

- Reduce the radiation dose and radiation emission time to the amount absolutely required. Any type of radiation has a detrimental effect on biological tissue.
- If an XRS 4 is mounted in the adapter plate of the stand and an applicator is mounted onto the mounting shaft of the adapter plate, this activates the optical safety interlock, thus enabling the INTRABEAM 600 to generate X-ray radiation.
- Suitable radiation shielding must be used to protect patient tissue outside the target area and OR personnel during use.
- Only work with the system after it has been properly calibrated. For notes on calibration see Page 39.
- Never operate the system unattended.
- Never cover any ventilation openings.

After each operation of the device

- Shut down the computer properly before switching the device off.
- Use the power switch to switch the device off. The power switch must always be switched off whenever the device is not in use.

Use of the INTRABEAM 600

In order to ensure the safety of patients and operating personnel, the corresponding requirements regarding selection of personnel, safety, training, verification, sterilization and other prerequisites must be complied with.

All XRS 4 units are calibrated by ZEISS at the manufacturer's plant. ZEISS supplies each XRS 4 with a new calibration file, see Page 146. An XRS 4 cannot be operated in the INTRABEAM 600 without this associated calibration file.

Only qualified personnel thoroughly trained in the use of the INTRABEAM 600 may use the INTRABEAM 600. Observe the instructions provided in the section "Information for operators and users", Page 27.

Based on the dose rate and the penetration depth, the INTRABEAM 600 calculates the total exposure time. Manual input of the exposure time is not provided for.

- As a precondition for using the INTRABEAM 600, you must have the output of the XRS 4 calibrated by the manufacturer, Carl Zeiss Meditec AG, based on a known reference which is traceable to a national standard.
- Switch the XRS 4 on only to treat a patient or for inspection purposes, making sure that the XRS 4 is properly connected to a component for inspection.
- Insert the XRS 4 only in the support of the stand and mount an applicator only if the XRS 4 is also supposed to emit X-ray radiation. You must use a suitable radiation shielding to protect all of the patient's tissue which is not to be treated and the operating personnel.
- Handle the ionization chamber with care and do not touch the sensor area. The ionization chamber supplied with the INTRABEAM 600 is a sensitive and fragile device.
- In order to ensure electrical safety, always connect the INTRABEAM 600 to protective ground before using it.
- Do not remove any components from the INTRABEAM 600. Doing so can cause malfunctions.

- When handling contaminated components of the INTRABEAM 600, always wear protective clothing such as gloves, an apron or a surgical gown and eye protection.
- Only equipment approved by ZEISS may be connected to the interfaces of the INTRABEAM 600.
- You must rotate the PAICH all the way around ($\geq 360^\circ$) the part to be measured in order to measure the decentration error exactly.
- Observe the instructions provided in the section "Handling, transport and shipping", Page 278.
- Observe the instructions provided in the section "Requirements for device use", Page 93.

Data protection and information security

The operator, his IT representative and every user of the INTRABEAM 600 must ensure that the corresponding laws and data protection regulations of the respective country are observed. The required technical and organizational prerequisites for this must be created so that

- the data in the INTRABEAM 600 is protected against misuse, i.e. that any personal, protected data cannot be collected, processed or used without authorization.
- personal data is protected, i.e. that any individual details about the personal or factual conditions of a certain or identifiable individual person, i.e. all data which can be directly assigned to a person (employee, customer, supplier), e.g. marital status, type of employment, religion, income, etc. cannot be passed on if it directly relates to a person.
- the data is protected during data processing, i.e. the data is protected while being stored (collected, recorded or archived), while being transferred (forwarded to parties outside the company), while being modified (changed in content, including anonymization and aliasing), and while being deleted or blocked (labeling to restrict further processing or use).

- the data is protected even during other use (e.g. in case of internal transfer).
- Data sent to third parties, outside of the responsible body, are protected by the respective recipient.

Archiving of data

- The INTRABEAM 600 can only store data. It is not intended to be used for archiving data.
- Data can be deleted by any subsequent user!
- In the event of a repair, all data are deleted from the computer. I.e. any data that you have not previously archived will be lost.

Each user is responsible for archiving his own data.

Data carriers must be free from viruses.

The operator is responsible for network security and ensuring that the external data carriers (USB sticks, CDs and DVDs) used are free from viruses.

Virus infestation of the software represents a hazard for the patient due to possible overexposure to radiation and/or prolongation of operating and anesthetic times.

Safety devices of the INTRABEAM 600

Optical interlock system

An optical interlock system has been incorporated into the XRS 4 to detect the type of radiation shielding incorporated in the device in which the XRS 4 may be installed.

If the XRS 4 is put into its holder and a ZEISS applicator is connected, the optical interlock system closes and allows the generation of X-ray beams. When X-ray radiation is generated, an audible acoustic signal is emitted from the loudspeaker of the control unit and the external warning lamp lights up.

The status of the interlock system is always displayed in the "Treatment" window, see the "Interlock Status" display. The four possible interlock states are shown in the table below.

With the exception of the XRS 4 holder and the applicators, the production of X-ray radiation is permitted only when devices with adequate radiation shielding are attached to the XRS 4.

Interlock status	X-ray radiation generation	Device
None	Deactivated	None
Non-shielded	Deactivated	PAICH
Shielded	Permissible	PAICH with installed ionization chamber (IC) PDA
Treatment	Permissible	XRS 4 and applicator are properly mounted on stand

EMERGENCY STOP, control lamps and external interlock

- The CC 600 has an Emergency Stop switch.
- The CC 600 has a loudspeaker:
 - When all interlock connections have been closed, i.e. the INTRABEAM 600 is in the Ready status, long beeps are heard.
 - During X-radiation emission, short beeps are heard.
- The XRS 4 is equipped with an indicator light which illuminates when X-radiation is being generated.
- The CC 600 visually indicates the generation of X-ray radiation. The CC 600 has indicator lamps for the following:
 - Power On, color: white
 - X-ray On, color: orange
 - X-ray ready, color: green
 - Pause, color: blue

The indicator lamps are described on Page 70.

- There is an external interlock connection on the reverse side of the workplace which can be used both as an external interlock connection and as a remotely controlled display for X-ray production.

For explanations on the external interlock connection and the X-ray display, see Page 300.

The INTRABEAM 600 is designed to protect patients and users in the event of electrical problems by diverting currents over conductive external surfaces which are connected to the facility electrical system ground by the power supply cord.

A potential equalization connector is provided on the rear panel of the workplace. This enables other active devices to be connected to the same ground potential as the INTRABEAM 600.

The workplace is electrically fused to reduce possible hazards due to malfunctions of the electrical system. When changing fuses, you absolutely must make sure that the replacement fuses have the proper rating.

Calibrations

In order to use the INTRABEAM 600, the user must have the following components calibrated:

- XRS 4
- IC ionization chamber

ZEISS recommends annual calibration.

The PTW company recommends calibration of the UNIDOS E dosimeter every two years.

- Follow the further instructions provided in the section "Calibration regulations", Page 254.

Safety during tests and verification

If you use the XRS 4 in connection with the PAICH for regulating the probe, no X-ray radiation will be generated and the high voltage of the XRS 4 will remain switched off.

If you use the XRS 4 together with the PAICH, inserted in the IC ionization chamber in order to measure the emitted X-ray radiation, X-ray radiation will be generated.

With the shielding of the IC ionization chamber, the PAICH offers sufficient radiation protection when used as described in these Instructions for Use.

If you use the XRS 4 together with the PDA in order to check the isotropy, X-ray radiation will be generated.

The PDA offers sufficient radiation protection when used as described in these Instructions for Use.

Safety during treatment

Physicians should carefully consider the clinical aspects when selecting patients for treatment with the INTRABEAM 600, since patients who preoperatively exhibit pronounced edemas, bleeding or severe calcification in the treatment field may possibly be unsuitable for irradiation with the INTRABEAM 600.

X-ray radiation is generated only when the XRS 4 is activated.

The CC 600 controls and monitors the radiation output during treatment with the help of the IRM internal X-ray monitor. When the prescribed radiation dose has been reached, the internal X-ray monitor stops the treatment.

The IRM internal X-ray monitor detects the amount of radiation the patient is exposed to during treatment.

Before treatment begins, the software calculates the total exposure time based on the entered dose and the penetration depth.

If the IRM internal X-ray monitor does not switch off the treatment (e.g. because its counting rate has decreased), the treatment time then serves as an abort criterion, i.e. if the originally calculated total irradiation time is exceeded by more than 10%, the timer function cancels the treatment.

- Be sure to use a sterile applicator when inserting the X-ray probe of the XRS 4 in a patient.
- In order to ensure that the treating personnel is not exposed to an excessive radiation dose, the ambient radiation should be monitored with a suitable measuring device during treatment.
- Use a portable radiation shield to protect the operating personnel.
- Switch the X-ray radiation off before removing the X-ray probe with the mounted applicator from the patient.

Heat build-up during treatment


One unavoidable side effect of generating X-ray radiation is the simultaneous generation of thermal energy in the X-ray probe of the XRS 4 and the dissipation of thermal energy into the patient's surrounding tissue.

- Please note: The surface temperature of applicators may exceed 41°C during the course of treatment. For further information, please refer to the Instructions for Use for the applicator in question.

Electrics of the INTRABEAM 600

The workplace is equipped with various electrical connectors on its rear panel. Only devices authorized by ZEISS may be connected to these connectors.

The workplace is electrically fused to reduce possible hazards due to malfunctions of the electrical system and to ensure sustained protection against fire.

If a fuse needs to be replaced, it is essential that you replace it with a fuse which has the same ratings and type. Observe the data on the rear panel of the respective device.

Use only power cords which comply with the standard IEC 60601-1 and maintain proper grounding of the INTRABEAM 600.

The INTRABEAM 600 must always be connected to the protective ground of the room when it is in use in order to ensure electrical safety. Always use the power cord that is supplied with the workplace and ensure proper grounding of the power supply.

ZEISS recommends additionally grounding the workplace by connecting the potential equalization to the stand.

The maximum voltage supplied to the XRS 4 equals ± 15 V.

Any high voltage hazard is prevented by generating the high voltage in a well-insulated housing, i.e. the electronic system is surrounded by the metal case of the XRS 4. The metal case of the XRS 4 is connected to room ground via the XRS cable and the workplace power cord.

- Take appropriate precautions when using other grounded or ungrounded devices in conjunction with the INTRABEAM 600.

The workplace has a potential equalization connector (1).

If the workplace is connected to the facility power supply (in compliance with IEC 60601-1), then there is a potential equalization via the power cord and the protective ground conductor of the facility power supply.

You can use the potential equalization connector (1) to connect other active devices to the same ground potential as the INTRABEAM 600. Electrical potentials between the INTRABEAM 600 and other grounded devices in the vicinity can be minimized by connecting the potential equalization connectors of all other devices to the potential equalization connector of the control unit.

The workplace can be connected to the power grid via the rubber IEC connector with a power cord.

- Always run the power cord through the coil (2) to prevent unintentional loosening of the plug connection.

The electrical components of the INTRABEAM 600 are not defibrillation-proof. Therefore, remove the XRS 4 from the patient before performing defibrillation.

Precautions regarding EMC

The INTRABEAM 600 complies with the EMC requirements of IEC 60601-1-2.

The INTRABEAM 600 must not be used directly adjacent to other ME equipment and ME systems. If use in the immediate vicinity of other devices cannot be avoided, you must constantly monitor the INTRABEAM 600 during operation to check compliance with its intended use.

- When operating the INTRABEAM 600, please observe the following precautions regarding EMC:
 - Use only spare parts approved by ZEISS for this device.
 - Do not use any portable or mobile RF communication equipment in the vicinity of the device, as it cannot be excluded this will impair the functioning of the device.
 - Do not use any cell phones near the device. They represent a potential risk for the proper functioning of medical equipment.
 - Malfunctions may occur, depending on a variety of local factors. They cannot be predicted and are not by any means calculable. Strong magnetic fields can have an impact on the XRS 4.
 - Do not expose the XRS 4 to strong magnetic fields.
 - Use the INTRABEAM 600 in hospitals only if no active RF surgical units are located in the vicinity.
 - Do not operate the INTRABEAM 600 in RF shielded rooms that are used for magnetic resonance imaging, as high-intensity electromagnetic interferences occur there.
- Please note the EMC guidelines in the section "EMC - Electromagnetic Compatibility", Page 303ff.

Fire, anesthetics and oxygen

You should not use inflammable substances such as anesthetics and oxygen in the vicinity of activated electrical components of the INTRABEAM 600 (CC 600, dosimeter, computer and monitor).

You must not switch on any electrically operated components of the INTRABEAM 600 if they located are in the vicinity (within 25 cm) of flammable anesthetics or an atmosphere enriched with oxygen or nitrogen monoxide.

Cleaning, disinfecting, sterilizing

- In order to protect patients and operating personnel against infections, use only components which have been carefully cleaned and disinfected during treatment.
- Use only sterile applicators.
- Always cover the stand with a sterile disposable cover (INTRABEAM Drape).
- Do not recycle any disposable products.
- Use only the approved procedures described here to process the components.

For further details and explanations on the approved cleaning and disinfection procedures, please refer to Page 282 ff.

The sterilization procedures for sterilizable applicators are provided in the separate, respective Instructions for Use of the applicator concerned.

Clean and disinfect by following the respective care instructions for INTRABEAM 600 components and accessories:

- Workplace
- V-guide (workplace)
- Monitor (touchscreen monitor)
Please observe the special instructions for the touchscreen on Page 283.
- Keyboard (membrane keypad)
- Mouse
- XRS 4
- PDA
- PAICH
- XRS and PDA/PAICH cables
- Stand
- INTRABEAM Water Phantom
- Dosimeter*
- IC ionization chamber*

* If the UNIDOS E dosimeter and the IC ionization chamber need cleaning, please refer to the corresponding Instructions for Use supplied by PTW, Freiburg.

The following applicators are suitable for steam sterilization:

- Spherical applicators and set with transport and storage container
- Surface applicator set
- Flat applicator set

Sterile disposable products include:

- INTRABEAM Needle Applicator
- Radiation shield, flat 12x17"
- INTRABEAM Drape

Disposal of consumables (disposable components)

Used consumables, for example drapes and needle applicators, are infectious waste which you must dispose of properly.

- Follow the appropriate regulations in place in your clinic.
- Observe the corresponding national regulations.

Disposal of components of the INTRABEAM 600

ZEISS takes back unused components of the INTRABEAM 600 in a clean condition.

- In the event of disposal, please contact your local representative at the ZEISS sales branch in your country.

INTRABEAM 600 product identification

Observe the symbols and labels on the device!

If you notice that any of these labels is missing or has become illegible on your device, please contact us or one of our authorized representatives. We will supply you with replacements.

Warning and information labels on the INTRABEAM 600 Workplace


Front side of workplace

INTRABEAM
600

- 1 Logo of the manufacturer
- 2 Product name
- 3 Connector designation "Ion Chamber" (ionization chamber)

Rear side of workplace

4 Rating plate

The rating plate contains the following information:

Manufacturer (company name)	Carl Zeiss Meditec AG
Manufacturer	
Address of the manufacturer	Goeschwitzer Strasse 51-52 07745 Jena, Germany
Device name	INTRABEAM 600
Fulfilled standard	IEC 60601-2-8:2010
Rated voltage range (115V)	110V - 125V
Rated voltage range (230V)	220V - 240V
Connected load	300VA
Line frequency range	50Hz - 60Hz
Maximum total weight with full configuration	Max. 155 kg
"Observe disposal regulations" symbol	
Indication of origin	MADE IN GERMANY
Ionizing radiation symbol	
CSA label	
CE labeling with number of notifying body	

For further technical data of the INTRABEAM 600, refer to the section "Technical data" on Page 293.


For more information on the "Disposal regulations", refer to the section "Disposal" on Page 288.


5 Product and device identification


This label contains the following information:

Reference number	 6405
Serial number	 640210xxxx
Manufacturing date	 YYYY-MM-DD
UDI (Unique Device Identifier)	
Device identification (DI)	(01)0404953906459
Product identification (PI)	(11)YYMMDD(21)640500XXXX

For further order data of the INTRABEAM 600, refer to the section "Order data" on Page 290.


6 SIP label

The SIP label contains the following information:

- Manufacturer (company name)
- Your possibilities for contacting the device manufacturer, i.e., phone number, fax number and e-mail address of the local contact of the national Carl Zeiss sales organization.
- SIP number
A unique identification number assigned to your device.

7 "Observe transport position" warning label

"Observe Instructions for Use" indicating label:

Observe the Instructions for Use or accompanying documents.

Illustration displays the workplace in the transport position.

- Always follow these instructions prior to every transport:

Do not deposit any components on the workplace; stow everything away in the appropriate drawers.

- 1) Close all drawers and the folding compartment. Once all drawers and the folding compartment have been closed, they can be secured by locking them in this position.
- 2) Rotate the V-guide so that its flat, closed side points upward.
- 3) Position the two monitor support arms flush above each other. Lock the monitor support arms in this position.

This symbol indicates a potential hazard which could result if the transport position is not observed.

- Always set the workplace to the specified transport position prior to transport in order to prevent damage.

8 Fuses F1 and F2


This label on the power supply module informs you about the required device fuses. Change the line fuse of the workplace, see Page 277.

F1/F2 250V ~ T6,3 A-H

9 "Risk of crushing" warning label

Fingers may be crushed between the monitor and the two monitor support arms. Do not touch this area while the monitor is being moved.

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Workplace connector panel

10 Label: "Preset voltage range"

Specification of the line voltage range this device was preset to at the factory.

Symbol: "Caution": Check setting prior to initial operation.


11 Designation of connectors:

- LAN
- Ext Interlock
- "Potential equalization" symbol

12 "Observe the Instructions for Use" label

Observe instructions for use or accompanying documents.

13 "WARNING" label

The "WARNING" label contains the following information:

a) Warning text:

Grounding reliability can only be achieved when the equipment is connected to the receptacle marked "hospital only" or "hospital grade".

b) Graphic symbol: "Caution"

Observe the accompanying documents.



c) Ionizing radiation symbol

d) Warning text:

This X-ray unit may be dangerous to patient and operator unless safe exposure factors and the operating instructions are observed.

(In order to prevent the x-ray unit from posing a risk to patients or personnel, the guidelines on the safe administration of radiation and the Instructions for Use must be observed at all times.)

Workplace drawers and folding compartment

- 14** Label indicating maximum permissible load of drawer (sliding keyboard shelf): 2 kg

- 15** Label specifying contents of drawer: "Keyboard and mouse"
The PC keyboard and mouse can be operated on the drawer (sliding keyboard shelf) when it is open.

- 16** Label indicating maximum permissible load of drawer: 7 kg

- 17** Label specifying contents of drawer: "PDA/PAICH and XRS 4"
In addition to the PDA, PAICH and XRS 4, the ionization chamber also can be stored in this drawer.


- 18** Label indicating maximum permissible load of drawer: 7 kg

- 19** Label specifying contents of drawer: "Book"
The Instructions for Use are stored in this drawer.
A second XRS 4 also can be stored here if necessary.


- 20** Symbol indicating contents of folding compartment, i.e. "computer and CD/DVD":
The computer (with CD/DVD drive) and the dosimeter are stored in this folding compartment.

Inscriptions, labels and symbols on the components

Labeling of the XRS 4


Key to symbols

This graphic symbol is used to indicate the emission of X-ray radiation.

This graphic symbol is used on radiosurgical equipment to identify the X-ray source assembly.

"Calibrated" label: This label identifies X-ray sources that were calibrated with flat and surface applicators.


Identification of XRS 4 cable

- 1 Color coding: green
- 2 Classification icon
- 3 Engraving of connector manufacturer

Identification of PDA (photo diode array)

- 1 Rating label
- 2 Technical data for shielding effectiveness of PAICH

INTRABEAM

PDA Photo Diode Array

Part 304534-7000-003 Serial XXXXXX

1 —————

Sterilize before use


Carl Zeiss

Germany

2 ————— 2,0mm Pb Equivalent (50 KeV)

Identification of PDA cable

- 1 Color coding: yellow
- 2 Engraving of connector manufacturer

1

2

Identification of ionization chamber with shield (IC shield)**Inscription on ionization chamber**



- 1** Protective cap for sensor window (membrane) with inscription: "Remove protective cap before using"
- 2** Ionization chamber with inscription: "Caution: Do not touch membrane"

Inscription on shield

- 3** Technical data
- 4** Rating label for shield


Identification on connector

- 5** Certificate number of manufacturer "PTW"
- 6** Identification label of manufacturer "PTW"

Identification of the PAICH (probe adjuster and ionization chamber holder)

- 1 Rating label
- 2 Technical data for shielding effectiveness of PAICH

INTRABEAM

**PAICH Probe Adjuster
Ion Chamber Holder**

Part 304534-6000-058 Serial XXXXXX

1 —————

Sterilize before use


Carl Zeiss

Germany

2 ————— 2,0mm Pb Equivalent (50 KeV)

Label on battery compartment of UNIDOS E

This label means that no batteries or rechargeable batteries must be used, if UNIDOS E is applied together with INTRABEAM system.

System overview

System description	64
Overview	64
INTRABEAM Workplace	66
Connectors on the workplace	68
Principal components	69
Control Console 600 (CC 600)	69
Computer of the INTRABEAM 600	72
XRS 4 with IRM internal X-ray monitor	74
Protection components	77
Probe protector	77
V-guide	78
Verification devices	80
PAICH probe adjuster and ionization chamber holder	80
Ionization chamber "IC shield"	81
PDA (photo diode array)	83
UNIDOS E dosimeter	84
Accessories for your INTRABEAM 600	86
NC32 INTRABEAM Floor Stand	86
Applicators	88
INTRABEAM Water Phantom	89

System description

Overview

All components of the INTRABEAM 600 are combined in a single mobile workstation, the INTRABEAM Workplace.

With the "NC32 INTRABEAM Floor Stand" as a support for the "XRS 4" X-ray source, the INTRABEAM system is complete.

Special applicators are available for different treatments.

Main components of the INTRABEAM 600

The main components of the INTRABEAM 600 are integrated or stored in the workplace.

- INTRABEAM Control Console 600 (CC 600) for control of the X-ray source (XRS 4)
- Computer with keyboard, mouse and touchscreen monitor for configuration and operation of the INTRABEAM 600
- X-ray source (XRS 4) with IRM internal X-ray monitor for irradiation, positioned with the help of the stand

For illustrations, see Page 15.

Components for quality assurance

- PDA photo diode array (PDA)
 - Probe adjuster and ionization chamber holder (PAICH)
 - V-guide and XRS mount
 - Ionization chamber with shield (IC)
 - UNIDOS E dosimeter
- For illustrations, see Page 14.

Protection of components of INTRABEAM 600

- INTRABEAM Workplace
for safekeeping of the individual components and their safe transport inside the company.

- Probe protector
for X-ray probe of XRS 4; for illustration, see Page 77.

Accessories for the INTRABEAM 600

Since the products are not approved in all countries, please contact your local sales partner with regard to approvals.

For comprehensive information on the various accessory components of the INTRABEAM 600, refer to the respective separate Instructions for Use. For an overview, see Page 30.

- Stand (NC32 INTRABEAM Floor Stand)
- INTRABEAM Drape
- INTRABEAM Spherical Applicators
- INTRABEAM Needle Applicator
- INTRABEAM Flat Applicator Set
- INTRABEAM Surface Applicator Set

Accessory for testing (option)

INTRABEAM Water Phantom

INTRABEAM Workplace

The INTRABEAM Workplace contains all essential components of an INTRABEAM 600.

The workplace is used to store and internally transport an INTRABEAM 600 and to provide its power supply. The workplace serves as a system rack and control station for inspection of the INTRABEAM 600 and for the treatment process.

The workplace has been designed exclusively to store the components of an INTRABEAM 600.

The workplace can be moved manually on its four twin casters for internal transport. The wheels are equipped with cable guards. Two wheels are equipped with locking devices. The casters are ESD-sensitive.



The "NC32 INTRABEAM Floor Stand" is provided as an accessory for reliable positioning of the XRS 4.

To avoid the spreading of germs and infections, the workplace must not enter the sterile surgical area during treatment.

- 1 Touchscreen monitor
- 2 CC 600
- 3 Keyboard & mouse
- 4 Computer
- 5 Dosimeter
- 6 XRS cable
- 7 PDA/PAICH cable

Stored in the workplace:

- 8 Instructions for Use
- 9 PAICH
- 10 PDA
- 11 Ionization chamber with shielding
- 12 XRS 4

Connectors on the workplace

Several electrical connectors are located on the workplace.

On the front side of the workplace:

13 XRS cable, color coded "green"

14 PDA/PAICH cable, color coded "yellow"

15 Connector for ionization chamber


On the rear side of the workplace:

16 Power cord connector

17 Potential equalization connector

18 LAN (network connector)

19 External interlock

Principal components

Control Console 600 (CC 600)

Functionality of the CC 600

The CC 600 receives task packages from the computer for verification and during treatment.

The CC 600 implements these task packages and monitors the process independently. At the same time, it displays the current data.

The CC 600 supplies raw data to the computer, which uses this data to generate treatment and verification records.

The CC 600 comprises the following major components:

- Software and electronics to monitor and control various system components
- Emergency Stop button
- Display for setting parameters
- Control lamps and buttons
- Power supply unit that supplies the low voltage required by the XRS 4
- Air pressure sensor

Controls, displays, connections

Control panel of the CC 600

1 Emergency Off button

The CC 600 can be switched off, i.e. the emission of X-rays can be stopped immediately at any time, with this large, red button. This may be necessary if a critical situation occurs during treatment.

The X-ray emission also can be stopped by actuating the <Stop> button (10) in the "Treatment monitoring" window (see "Treatment monitoring" window Page 224).

When a treatment is canceled (Stop) the INTRABEAM software automatically saves the treatment parameters.

Following a treatment discontinuation, you must reverify the XRS 4, unless you continue treatment of the same patient within 90 minutes. In order to continue treating the same patient, you must actuate the "Continue with previous patient" button (checkbox) in the "Treatment prescription" window, see Page 191.

In order to switch the CC 600 back on, you first must unlock the Emergency Off button (by turning it and pulling it out).

2 POWER ON indicator light (color: White)

Illuminated when the power supply of the CC 600 is switched on.

3 Control panel display

visualizes various data required for monitoring the components of the INTRABEAM core system. The contents of the control panel display vary according to the current procedure, see Items 8 - 15.

4 X-RAYS ON indicator light (color: Orange).

as on the XRS X-ray source.

Illuminated when the XRS is generating X-ray radiation.

5 READY indicator light (color: Green)

Illuminated when all preparatory activities have been performed and the system is ready to generate X-ray radiation.

6 PAUSE indicator light (color: Blue)

Illuminated if the generation of the X-ray radiation is temporarily deactivated and the system is in the PAUSE mode.


7 <PAUSE> button

You can use the <PAUSE> button to temporarily interrupt the emission of X-ray radiation. The treatment parameters remain intact. The "PAUSE" indicator light lights up.



In order to continue the emission of X-ray radiation, you must go to the "Treatment monitoring" window, see Page 222, and actuate the <Continue> (4) button. On the CC 600, continuation of the treatment (for example by actuating the <PAUSE>) button again is not possible.

8 - 15 The information shown on the control panel display of the CC 600 varies according to the current procedure:

Control panel display during treatment:

Control panel display during quality assurance (QA):

Computer of the INTRABEAM 600

Configure your INTRABEAM 600 and operate it with the computer.

The computer (1) is integrated in the workplace. It is housed in the lower compartment of the workplace.

The computer features:

- a membrane keypad and mouse (2),
- a touchscreen monitor (3),
- a USB port (4) guided to the outside on the workplace,
- a DVD/CD-ROM drive (5).


Switching the computer and touchscreen on/off:

Switch the computer and touchscreen on and off in connection with the entire workplace. The picture on the touchscreen (brightness, contrast etc.) is non-adjustable.

For full operation of the INTRABEAM 600 (online mode), the computer in the workplace is connected to Control Console 600 (CC600) and the dosimeter via RS232 interfaces.

With the software installed on the computer, you can:

- import calibration data via the CD/DVD drive
- collect, activate, deactivate and reformat the data of the components used,
- check the system,
- enter treatment parameters,
- verify the set treatment parameters,
- perform treatments,
- check treatments,
- document and save verifications and treatments,
- export data via the USB stick.

XRS 4 with IRM internal X-ray monitor

The XRS 4 is a linear accelerator used to generate X-rays.

The XRS 4 generates a high voltage from the low DC voltage supplied by the CC 600, thus transmitting an electron beam to the X-ray probe. The high voltage generator in the XRS 4 is completely surrounded by a grounded metal housing so that injuries of the patient due to high voltage can be ruled out. Since the XRS 4 generates the high voltage itself, no high voltage cable leading to the XRS 4 is required. For operation, the XRS 4 is connected to the CC 600 with the XRS connecting cable (a low voltage cable).


The electron beam strikes a target in the probe tip, thus creating the X-ray radiation. The XRS 4 emits X-ray radiation with an almost spherical radiation field on the probe tip.

The probe tip is made of beryllium (Be), a material which is permeable to X-rays. The outer surface of the probe tip is coated with a protective layer.

You must have the XRS 4 calibrated regularly. Have the XRS 4 calibrated by ZEISS once a year. If national regulations require shorter calibration intervals, these must be observed.

The XRS 4 cannot be sterilized. You must not insert an X-ray probe directly into the patient's tissue. Always mount an applicator on the X-ray probe.

- 1 Cable connector of XRS 4
- 2 Shaft (deflector unit)
- 3 Probe tip (point of origin of X-ray radiation)
- 4 X-ray probe
- 5 X/Y marking to align PDA photo diode array
- 6 IRM internal X-ray monitor

IRM internal X-ray monitor

The IRM internal X-ray monitor is the main monitoring system for the radiation emitted to the patient.

The X-ray radiation generated at the tip of the X-ray probe has a spherical (ball-shaped) radiation pattern. A portion of the radiation thus returns to the XRS 4 via the X-ray probe. The IRM internal X-ray monitor detects this X-ray radiation in the XRS 4.

Before inserting the X-ray probe in the patient for treatment, you must check to see whether the IRM internal X-ray monitor is measuring a radiation rate. The radiation rate measured by the IRM internal X-ray monitor is then used during treatment as an indirect measure of the radiation emitted on the tip of the X-ray probe.

For more information on calibration of the X-ray power of the XRS 4, see Page 254.

10


Protection components

Probe protector

The probe protector (1) is used to protect the X-ray probe of the XRS 4 (4) against mechanical damage.

Whenever the XRS 4 (4) is located outside of its storage container and is not in use, you must slide the probe protector (1) over the X-ray probe by hand and fasten the locking screw (2) on the XRS X-ray source.

- 1 Probe protector
- 2 Probe protector locking screw
- 3 Applicator mounting shaft
- 4 XRS 4
- 5 Holder for XRS 4 (part of stand)

V-guide

The V-guide constitutes a combination of two guides. Area (1), the "V-guide" (principle of a V-groove guide), is used to mount the XRS 4, the PAICH and the PDA horizontally. Area (2), the XRS-mount, is designed according to the form of the XRS 4 housing and is used to mount the XRS 4 standing (upright). The V-guide is integrated in the working surface of the workplace. It is a removable unit. The V-guide can be inserted in the working surface of the workplace with its guide side facing down. The guides are then protected and the working surface on the workplace is larger.

V-guide, area (1)

The V-guide is used to mount an XRS 4 and a PAICH or a PDA such that the X-ray probe can be inserted in the inspection instruments safely.

Since the XRS 4 has a different housing size than the PAICH and the PDA, the V-guide is divided into two parts. Both areas of the V-guide are marked with symbols: (3) mounting area for an XRS 4 and (4) mounting area for a PAICH or a PDA.

The X-ray probe can be mechanically damaged or rendered useless when it is inserted in one of the inspection instruments. This can possibly lead to false irradiation, therefore:

- Always use the V-guide to insert the XRS 4 in a inspection instrument.
- Ensure that the XRS 4, the PAICH and the PDA are placed in the right area of the V-guide.


XRS mount, area (2)

The XRS mount is used to vertically mount an XRS 4 in order to safely set it up together with an applied inspection instrument in a vertical position (with the X-ray probe pointing upwards). The measurements are then performed on the XRS 4 in this vertical position.

To ensure that you can, for example, rotate the PAICH around the X-ray probe, you must place the XRS 4 with the PAICH in the XRS mount upright.

If you set up the XRS 4 without an XRS mount, it can fall over and thus be damaged or injure the user. Therefore:

- Always use the XRS mount to set up the XRS 4 in an upright position.

Verification devices

PAICH probe adjuster and ionization chamber holder

The PAICH probe adjuster and ionization chamber holder has the following main functions:


- Straightness measurement of X-ray probe
- In order to enable the user to center (straighten) the X-ray probe
- To provide a radiation-protected location for mounting an ionization chamber while the X-ray radiation of the XRS 4 is being measured.


In reference to the first two functions, the PAICH is sometimes referred to only as "probe adjuster" or "centering device".

Key

- 1 Light beam
- 2 LED
- 3 X-ray probe
- 4 Plunger
- 5 Photo detector

The illustration below shows a cross-section of the PAICH and X-ray probe.

The axis of the diode light source and the photodetector is positioned such that the edge of the X-ray probe partially interrupts the beam.

If the X-ray probe is straight, rotating the PAICH leads only to a minimal change in the light intensity that reaches the photodetector. If the X-ray probe is not straight, the intensity of the light signal changes when the PAICH is rotated.

The PAICH can be positioned such that you can straighten the X-ray probe with the pusher.

The following section on the ionization chamber describes how to use the PAICH as an ionization chamber holder.

1 Mount for ionization chamber

This mount contains a temperature sensor. The INTRABEAM 600 requires the exact temperature for the correct calculation of the dose during the "Output check", see "Dose rate output check", Page 182.

2 Pusher for straightening the X-ray probe

Ionization chamber "IC shield"

With the help of the ionization chamber, the dosimeter measures the dose rate of the XRS 4.

For quality assurance, you have to check the output of the XRS 4 and the IRM internal X-ray monitor prior to each treatment. In this way, the calibration data can be corrected and an exact dose delivery can be ensured.

The ionization chamber is a sensitive device and cannot be sterilized.

WARNING

Danger due to defective ionization chamber

If you touch the sensor window, you may damage the ionization chamber up to the point of malfunction. This could possibly lead to overexposure.

- Do not touch the sensor window, touch the ionization chamber only on its side or on its cable.

For proper operation and in order to ensure radiation safety, you must insert the ionization chamber in its IC shield correctly.


The ionization chamber supplied by ZEISS was calibrated by a certification body in units of gray per coulomb (Gy/C). This calibration complies with the requirements of the Physikalisch Technische Bundesanstalt (abbrev. PTB, National Metrology Institute of the FRG) and optionally with those of the US National Institute for Standards and Technology (NIST). The operator is responsible for ensuring that the ionization chamber is calibrated according to the applicable national regulations. If no national regulations apply, the operator should use a standard recognized in his/its country.

**WARNING****Danger due to lack of calibration**

If the calibration intervals of the IC ionization chamber are not observed, the quality assurance measurements may be faulty. This could possibly lead to under- or overexposure.

- Have your ionization chambers calibrated every two years as recommended by the PTW company.

- 1 Sensor window of ionization chamber
- 2 IC shield
- 3 Connecting cable, connection to dosimeter

PDA (photo diode array)

The PDA (photo diode array) is used to measure the isotropy of the radiation emitted by the probe tip of the XRS 4.

The PDA comprises 5 photo diodes which are centered on the four side surfaces and on the top surface of a cube such that each of them is equidistant from the radiation center point of the XRS 4's probe tip.

During quality assurance, the signals of the PDA are displayed in the "Dynamic offsets" window and used to measure the distribution (isotropy) of the X-ray radiation emitted by the probe tip. Based on these measurements of the XRS 4, the isotropy of the radiation field is set automatically.

1 Connection (cable to CC 600)

2 Opening for inserting the X-ray probe

UNIDOS E dosimeter

An INTRABEAM 600 comes standard with a dosimeter, model UNIDOS E, from the PTW company. The dosimeter is used to measure the current generated by the IC ionization chamber while the latter is exposed to X-ray radiation. The UNIDOS E dosimeter is suitable for continuous operation.


When the XRS X-ray source is emitting X-ray radiation, the IC ionization chamber typically generates currents of 1 to 60 picoamperes [pA] if a bias voltage of +300 V DC is applied. (This bias voltage is automatically set by the INTRABEAM 600.)

When using the UNIDOS E dosimeter from PTW, refer to its Instructions for Use. Please read through the Instructions for Use for the UNIDOS E dosimeter from PTW completely before using it.

The UNIDOS E dosimeter requires a warm-up phase of at least 5 minutes before you can begin with the measurements.

Accessories for your INTRABEAM 600

NC32 INTRABEAM Floor Stand


The "NC32 INTRABEAM Floor Stand" is a required accessory for your INTRABEAM 600.

This stand is required to position the XRS 4 with the mounted applicator during treatment. It can be balanced and safely positioned using magnetic brakes.

The stand is intended for use in hospitals, clinics or other human medical institutions.

Detailed information on this stand is available in the separate Instructions for Use of the "NC32 INTRABEAM Floor Stand", see list on Page 30. For order data, refer to list on Page 291.

- 1 Axis of motion
- 2 Connector panel
- 3 Stand base with locking devices
- 4 Control panel
- 5 Storage tray
- 6 XRS 4 (component of the INTRABEAM 600)
- 7 Handgrip with release button for magnetic brakes
- 8 Head of INTRABEAM Floor Stand

Applicators

INTRABEAM Spherical Applicators

The INTRABEAM Spherical Applicators are accessories of the INTRABEAM 600.

They are used to apply a prescribed radiation dose into the tumor bed or on tumor margins during intracavitary or intraoperative radiotherapy with the INTRABEAM 600.

Detailed information on these accessories is available in the separate Instructions for Use of the "INTRABEAM Spherical Applicators", see list on Page 30.

INTRABEAM Needle Applicator

"INTRABEAM Needle Applicator" Sets (comprising a needle applicator and a guide tube) are accessories of the INTRABEAM 600.

They are used in connection with the INTRABEAM 600 for the intraoperative irradiation of tissue, including the irradiation of intracranial tumors. Due to its particularly slim construction, tumors can be reached even in narrow cavities.

Detailed information on these accessories is available in the separate Instructions for Use of the "INTRABEAM Needle Applicator", see list on Page 30.

INTRABEAM Surface Applicator Set

The INTRABEAM Surface Applicator Sets are accessories of the INTRABEAM 600. They are only intended for delivering a prescribed radiation dose in applications with the INTRABEAM 600

- during intraoperative radiotherapy, on a surgically exposed surface or in a tumor bed and
- in the treatment of tumors on the body surface.

The surface applicator is designed such that a flat radiation field is generated directly below a circular application area.

Detailed information on these accessories is available in the separate Instructions for Use of the "INTRABEAM Surface Applicator Set", see list on Page 30.

INTRABEAM Flat Applicator Set

The INTRABEAM Flat Applicator Sets are accessories of the INTRABEAM 600. They are only intended for delivering a prescribed radiation dose in applica-

tions with the INTRABEAM 600

- during intraoperative radiotherapy, on a surgically exposed surface or in a tumor bed and
- in the treatment of tumors on the body surface.

The flat applicator is designed such that a flat radiation field is generated in water at a distance of 5 mm below its circular application area.


Detailed information on these accessories is available in the separate Instructions for Use of the "INTRABEAM Flat Applicator Set", see list on Page 30.

INTRABEAM Water Phantom

The INTRABEAM Water Phantom is an optional accessory used for the inspection of your INTRABEAM 600.

The water phantom is used to measure the depth dose curve and the isotropy of the XRS 4 to check whether the measured values are within the expected range.

Detailed information on this optional accessory is available in the separate Instructions for Use of the "INTRABEAM Water Phantom", see list on Page 30. For order data, see list on Page 291.

This image shows a single sheet of white paper with horizontal blue ruling lines. The lines are evenly spaced and run across the width of the page. There is no text or other markings on the paper.

Preparations for use

Requirements for system operation	93
Preparing the workplace for use	94
Transporting the workplace	94
Connecting the workplace to the power supply	97
Visual system check	97
Spatial arrangement of the treatment facility	99
Connecting	102
Workflow	104
Structure of the software user interface	106
Tab bar	107
Menu bar	108
Buttons	108
Startup of the INTRABEAM 600	111
Switching the workplace on	111
Starting the computer	112
INTRABEAM Start Manager (start window)	113
1 <INTRABEAM> button	114
Select working status	116
2 "radiance" software	117
3 "LINAC Configuration tool" software	118
4 Importing calibration data	119
5 Data Management	120
Menus in the menu bar	124
1 <File> menu	124
2 <Tools> menu	125

3 <QA Reports> menu	133
4 Menu item <Patient Database>	133
5 Menu item <About>	138
Equipment inventory	140
Equipment inventory, "XRS linear accelerator" group	142
Recalibrating the XRS 4	146
Equipment inventory, "Equipment" group	148
Equipment inventory, "Applicators" group	152
Quality assurance	158
Performing quality assurance	160
Video instruction, video window	166
XRS tip adjustment	168
Dynamic Offsets (adjusting the beam deflection)	173
Isotropy check	176
Output check of dose rate	182
Treatment prescription	190
Selecting a patient	190
Marking the treatment site	192
Entering treatment parameters	194
Planning fractional treatment	202
Treatment approval	205

Requirements for system operation

WARNING

System set for medical electrical systems

Any additional equipment connected to medical electrical devices must demonstrably comply with the applicable standards and directives (for example IEC 60950-1, modified + Cor. + A1, modified + A1/Cor. + A2, modified for devices that process data).

Furthermore, all configurations must comply with the normative requirements for medical electrical systems (see IEC 60601-1-1 or IEC 60601-1 + A1, section 16).

Anyone who connects additional equipment to medical electrical devices is a system configurer and as such is responsible for the system's compliance with the normative requirements for medical electrical systems.

Please note that local legislation takes priority over the above-mentioned normative requirements. If you have any questions, please contact your local dealer or ZEISS Service.

Space requirements

Adequate space for operation of the INTRABEAM 600 must be provided both outside the patient care area and in the treatment room during treatment.

The INTRABEAM 600 generates X-ray radiation and contains sensitive components. Vibrations and shocks can cause malfunctions of and damage to the sensitive components of the INTRABEAM 600.

- Make sure that the components of the INTRABEAM 600 always have a vibration-free position.

Strong magnetic fields can have an impact on the XRS 4:

- Do not expose the XRS 4 to strong magnetic fields.
- The INTRABEAM 600 is only to be used in hospitals if there are no active RF surgical units in the vicinity.
- Do not operate the INTRABEAM 600 in RF shielded rooms that are used for magnetic resonance imaging, as high-intensity electromagnetic interferences occur there.

The specified ambient requirements, see Page 311, must be met.

Electrical requirements

The line voltage required and the maximum power consumption to be expected for the INTRABEAM 600 (without stand) are specified on the rating plates and in the chapter on "System data".

The wall outlet must be equipped with a properly connected protective ground contact. Do not use any multiple sockets, but rather, only single sockets or an isolating transformer.

Preparing the workplace for use

Transporting the workplace

Transporting the workplace over long distances


- Transport the workplace over long distances, e.g. when relocating it or returning it for repair, only in its original packaging or in special return packaging. Please contact your dealer or ZEISS Service for this purpose.

Transporting the workplace inside the building

The workplace is used exclusively for transporting and storing the components of an INTRABEAM 600. The workplace serves as a system rack and control station for inspection of the INTRABEAM 600 and for the treatment process.

Prior to transport inside the building

- Heed the "Observe transport position" warning label on the rear panel of the workplace.

1 "Observe transport position" label

"Observe Instructions for Use" indicating label:

Observe the Instructions for Use or accompanying documents.

Illustration displays the workplace in the transport position.

- Always follow these instructions prior to every transport:

Do not leave any loose components deposited on the working surface of the workplace, in the V guide or in the XRS mount. Pack all components in the respective compartments of the drawers provided for this purpose.

- 1) Close all drawers and the folding compartment. Once all drawers and the folding compartment have been closed, they can be secured by locking them in this position.
- 2) Rotate the V-guide so that its flat, closed side points upward.
- 3) Position the two monitor support arms flush above each other. Lock the monitor support arms in this position.

This symbol indicates a potential hazard which could result if the transport position is not observed.

- Always set the workplace to the specified transport position prior to transport in order to prevent damage.

2 "Risk of crushing" warning label

Fingers may be crushed between the monitor and the two monitor support arms. Do not touch this area while the monitor is being moved.

- Disconnect all external cable connections properly and place the cables in the workplace holding fixtures provided for this purpose.
- Do not let any cables lie on the floor loosely.

CAUTION

If the workplace is loaded improperly or beyond its permissible limits, it may tip over, thus injuring the user.

- Load the workplace only with the components of an INTRABEAM 600.

The workplace was not designed to accommodate persons climbing on its drawers and shelves or sitting on it.

Transport inside the building

The workplace can be moved manually on four dual wheels. The wheels are equipped with cable guards (1). All wheels are equipped with locking devices (2).

- Release the locking devices of the wheels by raising the tabs.
- Push or pull the workplace manually only by grasping one of the side hand rails provided for this purpose.
- Avoid all jolting and do not move the workplace over steps or raised landings.

CAUTION

If the workplace is moved on inclined surfaces, it may tip over, thus injuring the user.

- Do not move the workplace over inclines which slope more than 5°.

Workplace in a new position

NOTE

The casters of the workplace will be damaged if they are left standing in a disinfection bath for a longer period of time.


- Wipe the disinfectant off of the wheels immediately.

CAUTION

If the workplace is placed in the immediate patient area (sterile operating area), germs and infections may be spread..

- Ensure during treatment that the workplace does not enter the immediate sterile patient area.
- Do not park the workplace on an incline.
- Secure the workplace with the locking devices (2) of the wheels by pressing down on the respective tabs.
- Make sure that the workplace is standing securely and horizontally and cannot roll away.
- Always lay all cables so that they do not obstruct the user.

Connecting the workplace to the power supply

- Operate the workplace using only the supplied power cord.
- Ensure that the workplace can be disconnected from the power grid quickly in an emergency. I.e. position the workplace at the location where it is used so that the power switch on the back of the workplace always remains accessible for operation or connect the power cord of the workplace to a building wall outlet that always remains freely accessible.

The line-powered components in the workplace (CC 600, computer and dosimeter) are run via the integrated isolating transformer. The workplace has a connector for potential equalization.

Visual system check

- Start preparing for device use by inspecting the individual components of your INTRABEAM 600.
- Ensure that all components are available completely and in sufficient quantity.

WARNING

Do not use any unchecked components!

If components of the INTRABEAM 600 are not checked prior to use (visual inspection), defective components may be used unnoticed, thus causing faulty measuring results which could possibly lead to false irradiation.

- Ensure that the individual components do not show any obvious signs of damage or modifications such as, for instance:
 - Deformation, dents, scratches
 - Contamination, discolorations
 - Kinked, broken or frayed cables

**WARNING****Use only components which have been approved by ZEISS!**

If the user uses unapproved components for measurement (instead of the IC ionization chamber and the UNIDOS E dosimeter), this can result in faulty measurements which could possibly lead to false irradiation.

- Ensure that only components which have been approved by ZEISS are used (visual inspection).
- Ensure that the components have been calibrated.

The use of scratched or otherwise damaged data carriers (CDs or DVDs) with calibration files can cause problems with data that cannot be imported.

- Ensure that the data carrier is not damaged before you use it.

NOTE**Malfunctions due to dirty interlock surfaces!**

Dirty interlock surfaces can lead to situations in which the interlock is not closed and the INTRABEAM 600 cannot generate X-ray radiation.

- Ensure that the interlock surfaces are clean before each use (visual inspection).

Spatial arrangement of the treatment facility

A proposal for the spatial arrangement of equipment and personnel during treatment with the INTRABEAM 600 is shown in the drawings on Page 100ff.

The treatment room should be divided into a sterile and a nonsterile area, as shown in the general plan.

Sterile area

The sterile area is a sterile field which offers room for the patient, the surgical team and those components of the INTRABEAM 600 which directly or indirectly come into contact with the patient.

If unsterile materials and/or components are brought into the sterile operating area, that can cause cross-contamination in this operating area.

- Ensure that only sterile materials and/or components are brought into the sterile operating area.

CAUTION

Failure to use an INTRABEAM Drape can cause cross-contamination in the sterile operating area.

- Always use the stand in the operating room only with an INTRABEAM Drape.

A person located in the sterile field must attend to the handling and operation of the sterile components of the INTRABEAM 600 while setting up the equipment and during the verification processes before treatment.

This person must be familiar with the required tasks for inspection of the components of the INTRABEAM 600 for correct functioning and be able to carry out these tasks properly using aseptic techniques.


Arrangement for a tumor bed treatment,
e.g. for breast cancer

- 1** OR, sterile area
- 2** OR, nonsterile area
- 3** Patient
- 4** XRS 4 with sterile applicator
- 5** Stand (with INTRABEAM Drape)
- 6** Surgical assistant
- 7** Surgeon
- 8** Anesthetist
- 9** Surgical instruments
- 10** OR nurse
- 11** Radiotherapy personnel (physicists),
operation of the CC 600 and the computer

- 12** Workplace,
always to be positioned outside of the sterial area.

The workplace contains the complete equipment of an INTRABEAM 600:

- PDA
- PAICH
- INTRABEAM Drape
- Cables
- CC 600
- Monitor
- Computer
- Dosimeter and ionization chamber
- V-guide

Connecting

NOTE**Danger due to defective connectors or cables**

The use of defective connectors or cables can lead to malfunction or damage of the INTRABEAM 600.

- Before using the INTRABEAM 600, perform a visual inspection of all plug connections and cables for possible damage.
- If you detect any damage, arrange for it to be repaired by ZEISS Service.

Incorrect connections can lead to function failure or damage of the INTRABEAM Core System.

- Never attempt to forcefully connect any electrical connectors (plugs, sockets). If a connector does not fit a socket, check if they are compatible with each other. The male and female components are mechanically coded.
- Do not pull at the cables.
- The plug connections are self-locking. Slide back the respective sleeve to separate the connections.

NOTE**Danger due to connection of unauthorized devices**

If unauthorized devices, for example USB sticks, hard disk drives or printers, are connected to the USB port, this may lead to a situation where the driver software of the unauthorized devices installs itself, thus interfering with the INTRABEAM software.

- Connect only USB sticks which are approved by ZEISS.

NOTE**Virus-free!**


The user is responsible for ensuring that the USB medium used for data exchange is virus-free.

NOTE**Secure network!**

The user (IT manager) is responsible for ensuring that the network is connected to the system via a suitable security design (e.g. firewall) to protect it against unauthorized access from outside.

Connecting the workplace

- Make the following cable connections:
 - Workplace to power grid
Run the power cord through the coil to prevent unintentional loosening of the plug connection.
 - Establish potential equalization, connection for potential equalization
 - Workplace with XRS 4
XRS cable, color coded "green" to XRS 4
 - Workplace (CC 600 control unit) with PAICH or PDA
PDA/PAICH cable, color coded "yellow", to PDA or PAICH
 - Workplace with ionization chamber
Cable of ionization chamber to "Ion Chamber" connector of workplace
 - Workplace with data network, LAN connection (Ethernet port) for network connection of your INTRABEAM 600. The requirements for connecting the device to an existing IT network (according to IEC 60601-1) can be found on Page 298.
 - Connect external interlock (option)


Workflow

If the steps required for quality assurance are performed improperly, this can cause damage to the X-ray probe and also result in the INTRABEAM 600 switching off and not allowing irradiation with the respective XRS 4.

- Always adhere to the prescribed workflow when working with the INTRABEAM 600.

Starting up the INTRABEAM 600

1. Switch on workplace
2. Log into computer
3. INTRABEAM Start Manager
4. Start INTRABEAM software

Structure of the software user interface

The user interface is divided into multiple task-related windows. For details, diagrams and reports, you can change to the window of a different task area or open a trailing window.

All control elements (buttons) on the touchscreen can be pressed or clicked with the mouse.

You can enter numerical values and texts (e.g. comments) using the PC keyboard or via the virtual keyboard on the touchscreen.

Key functions in the header areas of the windows

In the header areas of windows you will find:


1 Menu bar

The individual menus of the menu bar have a pushbutton function. If you actuate one of these menus (the text), the corresponding menu opens to offer you the individual menu items (lines of the list) for your selection. The individual menu items also are buttons. If you actuate a menu item, the corresponding window is displayed or the corresponding function is executed.

2 Tab bar

The individual tabs (resembling index card tabs) have a pushbutton function. Actuating a tab opens the corresponding window.

Tab bar


The tab bar contains all functional groups of the software in their logical order, arranged from left to right. I.e. in a regular workflow, you call these functional groups in the given order.

When you click on a tab, you go to the window of the corresponding functional group:

- Equipment inventory
Collect data of existing equipment, see Page 140.
- Quality assurance
Inspection of X-ray source, see Page 158.
- Treatment prescription
Entry of treatment parameters, see Page 190.
- Treatment approval
Inspection of treatment parameters, see Page 205.
- Treatment monitoring
Performance of treatment, see Page 244.
- Treatment review
Display of treatment parameters, see Page 231.
- Calibration
Calibrations, see Page 257.

Menu bar

The menu bar contains the following menus:


- 1 <File>
For more information, see Page 124.
- 2 <Tools>
For more information, see Page 125.
- 3 <QA Reports>
For more information, see Page 133.
- 4 <Patient Database>
For more information, see Page 109.
- 5 <About>
For more information, see Page 138.


Buttons

Function buttons

The underlying function is triggered by actuating the button.

Since the control elements of the software user interface can be operated either with the mouse or via keyboard in the INTRABEAM 600, instructions for operation procedures are referred to in a neutral manner with "actuate" in the Instructions for Use.

Keys with status display

The underlying function is triggered by actuating the button.

Color code for status display

- Yellow** The task concerned still must be completed or a certain component e.g. is ready for verification.
- Green** The task concerned has been properly completed or a certain component e.g. has been successfully verified.
- Red** Error! The task concerned was not properly completed or a certain component cannot be used, e.g. because it did not pass the verification.

Menu buttons

Menu buttons are marked by a triangle symbol. If you actuate a menu button, the corresponding menu opens. In the opened menu, you can actuate the individual menu lines, so that the relevant window opens or the relevant function is executed.

Input fields

*Patient ID	*Family name	*Given name
01294870	Maier	Hanz

You can enter texts and numerical values in input fields either via the PC keyboard or via the keypad on the touchscreen.

Input fields are titled to inform you concerning the type of input (text and/or numerical value) required.

If the title of the input field is marked with an * (asterisk), it is a mandatory field. I.e. an entry in this field is essential for proper execution of the corresponding software function.


The following input rules apply to all inputs of texts and numerical values:

- Do not use any special characters such as < > | * ? : / \
- Do not enter more than 20 characters in an input field.

Input fields for search function

A search can easily be performed with the help of the wildcard function, i.e. partial entries of texts or numbers preceded or followed by an * (asterisk).

Display and comment fields

In some windows, different types of fields can be displayed at the same location. Using the keys above the fields, you can switch back and forth between the different types.

<Log File>

This is a pure display field, i.e. the current process is automatically logged. You cannot edit anything here.

<Comments>

You can enter notices and comments in this field via the keyboard.


Scroll function

A scroll bar is displayed on the right next to some fields (1). If the contents of the field are extensive, not all of them can be displayed at once. However, you can view the entire contents of the field by scrolling them up and down with the scroll bar.

Startup of the INTRABEAM 600

The line-powered components integrated in the workplace (computer, CC 600 and dosimeter) can be connected to the power grid only via the isolating transformer of the workplace.

Switching the workplace on

Using the power switch (1) on the back of the workplace, you can simultaneously switch all line-powered components of your INTRABEAM 600 (computer, CC 600 and dosimeter) on and off.

In position "I" the workplace is connected to line voltage.
In position "O" the workplace is separated from the line voltage.

Briefly switch the workplace on and off.

After the workplace is switched on, the computer initiates its start routine: If the workplace is switched off again immediately, the start routine of the computer is canceled. Switching the workplace on again immediately can damage the computer. Observe the following information:

NOTE

Risk of damage to the computer

Briefly switching the workplace on and off can damage the computer.

- After you shut down the workplace, wait for at least 30 sec. before switching it on again.

Starting the computer

Once you have switched the workplace on, the computer initiates the start routine:

- The computer login window is displayed.
- To start the computer software, You must enter the user ID and the corresponding password.

You require a user ID and the corresponding password to use the computer. The passwords for the computer login and the administrator login are assigned by ZEISS Service.

The administrator sets up the access authorization with password for the user.

While the computer startup routine is running,

- the CC 600 performs a self-test
- and the dosimeter starts a warm-up phase. The dosimeter requires a warm-up phase of at least 5 minutes before you can begin with the measurements.

Computer software login

- User name
Enter your user ID for the computer.
The user ID for the computer is assigned by ZEISS Service.
- Password
Enter the password for the computer.
The password for the computer is assigned by ZEISS Service.
- Actuate the key with the arrow symbol or the <Return> key on your PC keyboard to execute the computer login.

If the login process is canceled and the computer switches off, observe the following:

NOTE

Risk of damage to the computer

- After the computer shuts down, wait at least 30 sec. before restarting it.


Once the computer software has started, the INTRABEAM Start Manager (start window) is displayed.

INTRABEAM Start Manager (start window)

As soon as you have executed the computer login using the login window, the "INTRABEAM Start Manager" window appears.

- Actuate the relevant button in order to go to one of the following functions:

- 1 <INTRABEAM>
For more information, see Page 140.
- 2 <radiance>
For more information, see Page 117.
- 3 <Linac Configuration>
For more information, see Page 117.
- 4 <Calibration Data Import>
For more information, see Page 119.
- 5 <Data Management>
For more information, see Page 120.


Return to "INTRABEAM Start Manager":

- Actuate the "Home" symbol (6) in the task bar.

1 <INTRABEAM> button

Start of INTRABEAM software:

- Actuate the <INTRABEAM> button. The "Software Login" window is displayed.

- 1 "User ID" input field**
Enter your user name (user ID).
The user ID is assigned by the user-side administrator.
- 2 "Password" input field**
Enter your password for the INTRABEAM software.
The initial password is assigned by the user-side administrator.
- 3 <Working status> menu button**
Select your desired working status:
Actuate the <Working status> menu button. The menu items <Offline> and <Online> are then offered to you for selection. For more information on the working status, see Page 116.
- 4 <Log in> button**
Once the user ID and password have been entered correctly and a working status has been selected, the INTRABEAM software can be started by actuating this button.
- 5 <Cancel> button**
If you actuate the <Cancel> button, all of your current entries and your selection in the <Working status> menu will be cleared.
- 6 <X> button**
If you actuate the <X> button, this stops the login process and returns you to the INTRABEAM Start Manager.

User name and password

During the initial installation of an INTRABEAM 600, ZEISS Service and the operator jointly define the basic settings for the name and password regulation.

In this process your username is used for logfile entries.

User groups

Five user groups are designated in the software:

Administrator	Physicist	Physician	Therapist	Service
---------------	-----------	-----------	-----------	---------

During the verification of the treatment parameters, the physician and physicist user groups must confirm these values. Ensure that only medical physicists are assigned to the physicist user group and only radiation oncologists are assigned to the physician user group.

Administrator

The administrator can only administer users. He is not able to operate the software.

During the initial installation of an INTRABEAM 600, ZEISS Service and the operator jointly define the name and password regulation for the administrator.

User

Within his own user group, any user can create and administer other users with the same access rights.

When the software is started, the window for user administration ("Administration" window) can be opened via the following path:
Menu bar, <Tools> / <Administration>.

The "Administration" window opens and the following functions are available for your selection in a pull-down menu:

New Account	Overwrite password	Delete access
Change password	List Users	Preferences

For more information, see Page 125.

Select working status

When you have entered your user name (1) and your password (2) during the software login, you then must select your desired working status:

- Actuate the "Working status" menu button (3); the menu items <Offline> and <Online> are then offered to you for selection. More information on the working status can be found below.

Offline status, <Offline> button (1)

"Offline" means working only with the computer. Only administrative work can be performed in this status. You can:

- administer the components of your INTRABEAM 600 (your inventory), i.e. work with the "Equipment inventory" window,
- view and administer patient information, i.e. work with the patient database,
- plan treatments, i.e. work with the "Treatment prescription" window,
- administer treatment data, i.e. work with the "Treatment review" window and the QA records.

Online status, <Online> button (2)

The online status is the full status. When the online status starts, the CC 600 and the computer are synchronized in that the CC 600 imports the current time and date from the computer.

The complete software can be used in the online status, i.e. you can:

- Administer the components of your INTRABEAM 600, i.e. work in the "Equipment inventory" functional group. To do this, actuate the corresponding tab.
- Perform quality assurance, i.e. work in the "Quality assurance" functional group in order to check the XRS 4 currently planned for use in a treatment. To do this, actuate the corresponding tab.
- Enter patient information and treatment parameters. Enter the patient information and treatment parameters for a pending treatment, i.e. work in the "Treatment prescription" functional group. To do this, actuate the corresponding tab.

- Have the treatment parameters entered checked by the radiotherapist (physician), i.e. work in the "Treatment approval" functional group.
- Perform and monitor the treatment, i.e. work in the "Treatment monitoring" functional group.
- Administer treatment data, i.e. work in the "Treatment review" functional group.
- In the "Calibration" functional group you can:
 - Determine the sensitivity of the beam deflection
 - Determine the linearity of the IRM internal X-ray monitor
 - Calibrate the air pressure sensor
 - Calibrate the temperature sensor

2 "radiance" software

radiance "radiance" (GMV Spain) is a software which is used for irradiation planning and simulating treatments performed with intraoperative radiation therapy devices. The version of radiance which, depending on the country-specific approval, is supplied with the INTRABEAM 600, supports 3D simulation of irradiation for INTRABEAM applicators.

"radiance" allows the import of medical 3D image data, the segmentation of clinically relevant target volumes, irradiation simulation and dose calculation for INTRABEAM applicators and the compilation of reports.

You can tell whether radiance is available on your INTRABEAM 600 based on the appearance of your INTRABEAM Start Manager.

Function call

- Actuate the <radiance button in the INTRABEAM Start Manager; the "radiance Login" window is then displayed.

3 "LINAC Configuration tool" software

LINAC Configuration tool Using the "LINAC Configuration tool" which belongs to the "radiance" software, you can create models of INTRABEAM X-ray sources (XRS 4) and INTRABEAM applicators based on calibration files.

For more information on the use and operation of radiance and the LINAC Config tool, please refer to the "radiance" Instructions for Use included in the scope of supply.

NOTE

Risk of operating errors

- Before using radiance and the LINAC Config tool, it is essential for you to familiarize yourself with the Instructions for Use of this product.

You can tell whether "radiance" software is available on your INTRABEAM 600 based on the appearance of your INTRABEAM Start Manager.

Function call

- Actuate the <Linac Configuration button in the INTRABEAM Start Manager; the "Linac Configuration Login" is then displayed.

4 Importing calibration data

Calibration data can be copied from a CD or DVD to the hard disk drive of the computer.

NOTE

Risk due to hardware or data carrier defect

If the computer hardware or the data carrier (CD/DVD) is damaged, no new calibration data can be imported. This may possibly make it necessary to postpone a treatment.

- Always import calibration data prior to a treatment in the context of equipment inventory or by the beginning of quality assurance at the latest.

NOTE

Risk due to vibrations and jolts

Vibrations and jolts impacting on the computer while data are being copied to the hard disk drive may cause incorrect or incomplete data to be imported.

- Ensure that the workplace always has a vibration-free location.

This calibration data import is possible only when the INTRABEAM software is not currently being used, i.e. when no user is currently logged in. This new calibration data will initially be accepted by the INTRABEAM software only after it is rebooted, e.g. for quality assurance.

- Shut down the INTRABEAM software before importing the new calibration data. To do this, go to the <File> menu and actuate the menu item <Exit>. Alternatively, you also can exit the INTRABEAM software by entering the key combination <Ctrl>+<Q> on the PC keyboard.

The INTRABEAM Start Manager window is displayed.

- Insert a CD or a DVD with calibration data in the CD drive of the computer.
- Actuate the <Calibration Data Import> button.

The INTRABEAM software automatically imports the calibration data. Finally, an information window shows whether or not the calibration data import was successful.

The new calibration data are displayed for checking and verification during the next start of the INTRABEAM software.

5 Data Management

Using the "File browser" window you can access a selection of released directories of your INTRABEAM 600.

Function call

- In the INTRABEAM Start Manager, actuate the <Data Management> button; the "Data Management" window is then displayed.

Different directories

The following directories can be opened by actuating the corresponding button:

<Treatment reports>

Text documents in PDF format

<QA Reports>

Quality assurance reports and verification reports in PDF format

<Session log files>

Troubleshooting files for ZEISS Service

<Error log files>

Troubleshooting files for ZEISS Service

<Screenshots>

You can store your images here.

<Miscellaneous>

Currently still without function

<Database>

Currently still without function

<User management>

File containing all users of your INTRABEAM 600, file for data backup in BIN format

Selecting all folders:


- Actuate the button (checkbox) next to "Select all folders".

Selecting all files of a folder:

- Actuate the button (checkbox) next to "Select all files".

or

- Actuate the checkbox of the relevant button when opening an individual folder.

Editing folders

In the individual folders you can:

- Select individual files or all files.
- Export files to a USB stick.
- Delete files.

Selecting a file:

- Actuate a file icon; the relevant file icon is then highlighted by a frame.

Selecting all files:

- Activate the button (checkbox) next to "Select all files". All file icons are then highlighted by a frame.

Deleting files:


- Select the files you want to delete and actuate the <Delete> button; the selected files are deleted from the computer hard disk.

Exporting files:


Individual files or entire folders can be exported (copied) to a USB stick.

The USB port is located in the keyboard drawer, to the left of the PC keyboard.

- Connect a suitable USB stick to the USB port; the type of USB stick connected is indicated on the "Destination" menu button.
- Select the files you want to copy and actuate the <Export> button; the selected files are copied to the USB stick.


Menus in the menu bar

In the menu bar the name (text) of the individual menus is linked to a corresponding button function. If you actuate one of these menus, the corresponding menu opens to offer you the individual menu items (lines of the list) for your selection.

The individual menu items also are buttons. If you tap on a menu item, either the corresponding window will be displayed or the corresponding function will be performed.

1 <File> menu


The <File> menu contains the menu item <Exit>


The INTRABEAM software can be shut down by actuating the menu item <Exit>.

Alternatively, you also can exit the INTRABEAM software by entering the key combination <Ctrl>+<Q> on the PC keyboard.

Before the INTRABEAM software is actually terminated, a confirmation prompt follows.


2 <Tools> menu

The <Tools> menu contains the following menu items:

- <Applicator management>
- <XRS Usage> / <View>
- <Administration>
- <Set Control Console Volume>
- <Login>

The individual functions are explained below.

<Tools> / <Applicator management> menu item

You can manually document how often an applicator has been sterilized with the <Applicator management> function.

In radiotherapy, often a selection of applicators of different sizes are provided to choose from during the procedure. All these applicators are then regarded as "contaminated" and must be resterilized.

Only the applicator actually used during treatment, i.e. the applicator selected during the "Treatment prescription", is thus recorded by the INTRABEAM 600 and its number of sterilizations is automatically incremented.

For the other applicators that were not used but nevertheless are to be regarded as "contaminated", you need to manually increase the sterilization count. It is advisable to document the number of sterilizations performed on these applicators during the treatment.

- Actuate the menu item <Applicator management>; the "Applicator management" window is displayed.

The "Applicator management" window is a tool for monitoring the processed applicators of your INTRABEAM 600. This auxiliary tool does not relieve you of the responsibility of correctly documenting the number of processings for all applicators.

WARNING

Permissible number of processings exceeded

Once an applicator's permissible number of processings has been exceeded, its cleanliness, sterility and biocompatibility no longer can be guaranteed. The functioning of the applicator no longer can be ensured.

Defective applicators may result in the emission of a radiation dose which does not correspond to the planned dose. The consequences are underexposure or radiation damage due to overexposure. If the patient or user comes into contact with an unsterile or no longer biocompatible applicator, the consequence of this is allergic reactions and infections.

- Ensure that all processing cycles are counted and that no applicator which has undergone more than the permissible number of processing cycles is used. For the maximum permissible number of processing cycles, see the table below.

Applicator type	Maximum permissible number of processing cycles
Spherical applicator	100
Surface applicator	50
Flat applicator	50

Selecting applicators

The "Applicator management" window contains two tables for applicators. The left table with the heading "Available" lists all of the active applicators currently available in your INTRABEAM 600. The right table with the heading "To sterilizer" is provided for the applicators which must be processed.

- Select the applicator from the "Available" table on the left which you regard as "contaminated" and move it to the "To sterilizer" table on the right:
 - Actuate the respective line in the "Available table".
 - Actuate the ">>" arrow button (pointing to the right).

Individual applicators (individual lines) also can be moved from one table to the other via drag and drop with the mouse.

You can use the "<<" button (pointing to the left) to correct your selection.

Incrementing the number of usages:

- Actuate the <Confirm> button once.
The number of usages of the applicators is then incremented by one; see the "Usage" column of the "To sterilizer" table on the right.

In the "Applicators" group of the "Equipment inventory" window (see Page 152) you can see that the number of usages displayed has been incremented; the number of usages is only displayed, i.e. cannot be manually incremented, there.

As soon as the maximum permissible number of processings is reached for one applicator, this applicator is set to "Inactive" (see "Equipment inventory, Applicators Group"); it therefore no longer can be used. (If you would like to continue using an applicator contrary to the recommended maximum number of processings, you can reactivate it. However, you must then confirm this with your password, see Page 152.)

Applicator usage management

Applicators

Available


Type	S/N	Size (cm)	Usage
Spherical	A63003	1.5	0
Spherical	B63001	2.0	0
Spherical	B63002	2.0	0
Spherical	H63002	5.0	0
Flat	101021	1.0	0
Flat	102001	1.0	0
Flat	310020	2.0	0
Flat	151034	6.0	0
Flat	151035	6.0	0
Flat	152001	6.0	0
Flat	152002	6.0	0

To Sterilize

Type	S/N	Size (cm)	Usage
✓ Spherical	A61065	1.5	0
✓ Spherical	H63001	5.0	1

Please move any applicators to be sterilized from the left list to the right list. All applicators on the right will have their number of sterilizations incremented in inventory when 'Confirm' button is pressed. You may drag and drop single items or use arrow buttons to move single or multiple items.

Cancel Confirm


Menu item <Tools> / <XRS Usage> / <View>

You can display the usage data (hours of operation) for each XRS 4 of your INTRABEAM 600.


The complete irradiation time and the portion of that time that was treatment time are indicated in hours.

- Actuate the menu item <XRS usage/ <View>;
The "XRS Usage Summary" is then displayed in the "Report browser" window.
- To save the report, please only use the <Save PDF> button.

(Blank page, for your notes . . .)

[illegible]

Menu item <Tools> / <Administration>.

This function is used for user administration.

- Actuate the menu item <Administration>; the "Administration" window is then displayed.

Pressing the "Select" button then opens a menu from which you can call the various user administration functions.

New Account

You can add additional users to your user group.

Change Password

The user who is currently logged in can change his password.

Overwrite password

A function **only** for the administrator.

List Users

List of all users who are registered in the system.

Delete access

A function **only** for the administrator. Used to delete a user from the user list.

Preferences


Setting the signature rules:

If you remove the check mark from the input field, this means that the physician no longer has to sign the treatment data.



If you remove the tick from the input field, making it unnecessary for a second person (e.g. a doctor) to check the treatment data, you may be in contravention of the existing legal requirements.

DICOM also can be switched on or off here.


The camera function is intended for future upgrades.

Default Administrator


Menu item <Tools> / <Set Control Console Volume>


This menu item can be used to set the volume of the loudspeaker integrated in the CC 600.

- Actuate the menu item <Set Control Console Volume>; the window for setting the volume of the CC 600 is then displayed.

The volume setting will be reset to the factory setting when the software is restarted.


Menu bar, <Tools> / <Login>.

Renewed callup of the INTRABEAM software (see Page 140), for example in order to change users.

- Actuate the <Login> menu line or enter the key combination <Ctrl>+<L> on the PC keyboard.


3 <QA Reports> menu

The <QA Reports> menu contains the menu items:

- <QA Reports> / <Open>
- <QA Reports> / <Open Latest>


Menu item <QS Protokolle> / <Open>

Once you have performed the "quality assurance" routines (checks of the XRS 4 with PAICH and PDA) in your INTRABEAM 600, you can display all of these results, i.e. all of the QA reports stored in your software, as a list.


- Select the required QA protocol in order to open it for reading.

Menu item <QA Reports> / <Open Latest>

Once you have carried out the "quality assurance" routines, you can display the most recent results from the XRS 4 checks performed with PAICH and PDS, i.e. the QA report that was last saved in your software.

4 Menu item <Patient Database>

The menu item <Patient Database> is used for the administration of patient information and the associated diagnostic data. It also is possible to import patient information from DICOM.

- Actuate <Patient Database> in the menu bar.
The "Patient administration" window is then displayed.

The following input rules apply to all inputs of patient information:

- Do not use any special characters such as < > | " ? : / \
- Do not enter more than 20 characters in an input field.

"Patient administration" window

All patient information is displayed in a single table.

Search functions

- Actuate the <Filtern nach> menu button, all columns of the table are offered for selection as a menu bar.
- Select a table column.
- Enter your search criterion in the "Search for" field.

If the entry is incomplete, the search filter starts automatically. The only exception is the date of birth; a complete entry is required here.

- Actuate a line in the table to select the patient concerned.
The table line is color highlighted. The information of the selected patient is displayed in the fields below the table under "Patient information".

Editing patient information

Only the information of patients who have not yet been treated can be edited.

- Edit the patient information in the respective input field.
- Actuate the <Apply> button to import the changes to the selected patient information to the patient administration.

Add new patient

- Actuate the <Add new patient> button
to import the current entries to the patient administration as a new patient.

All fields marked with an * (asterisk) are mandatory fields, i.e. you are required to make an entry there.

- Fill in the "Patient information" fields (at least the required fields), then tap on the <Apply> button.

Information on the treatment site and diagnosis

In order to store information for the diagnosis of the respective patient, you must enter it in the "Diagnosis" input field.

A diagnosis must be created for each patient, otherwise the data for the patient in question cannot be saved.

Only the diagnosis can be entered in "Patient Administration". You can enter specific information on the treatment site in "Treatment Prescription", on the button "Mark site" under <Patient Information>..

- Use unique descriptions, as identification will be performed based on this "diagnosis" later on.
- Actuate the <Apply> button.

<Patient ID>
<Family name>
<Given name>
<Date of birth>
<Sex>
<Insurance ID>
<Other ID>
<Diagnosis>
<Status>

Add diagnosis

Several diagnoses can be stored for a single patient.

- Select the patient concerned.
- Fill in the "Diagnosis" field.
- Actuate the <Add diagnosis> button.
A new line is added to the "Diagnosis" column of the table for each additional diagnosis of the patient concerned.

If diagnoses have been added to a patient, the <Delete diagnosis> button also is activated, i.e. diagnoses of this patient can be deleted by pressing this button.

Deleting a patient's information

- Actuate the <Delete patient> button to remove the selected patient from the patient administration.

Exiting the patient database

- Actuate the <OK> button.
You then return to the window where you actuated the menu item <Patient Database>.

Alt patients

Filter by: Search for:

Show All

Patient ID	Family name	Given name	Date of birth	Sex	Insurance ID	Other ID	Diagnosis	Status
1	Mustermann	Alex	1970-01-01	Other			spinal metastasis	Unlocked
4	Mustermann	Daniela	1960-04-04	Female			non melanoma	Unlocked
1	Mustermann	Maria	1960-01-01	Female			mama carcinoma	Unlocked
2	Mustermann	Max	1965-02-02	Male			glioblastoma	Unlocked
5	Mustermann	Michael	1980-05-05	Male			H&N cancer	Unlocked

Patient information

Patient ID: 1 Family name: Mustermann Given name: Alex

Sex: Other Date of birth: 1970-01-01 Insurance ID: Other ID:

Site information: spinal metastasis

Buttons: Add new patient, Apply, Delete patient, Add diagnosis, Apply, Import patient, OK

Importing patient information from DICOM

You can copy personal information (without diagnostic information) to the patient database of your INTRABEAM 600 via the DICOM interface.

Sequence:

- Search: first you must search for a patient via DICOM.
- Select: select the required patients from the search results.
- Import: import the selected patient data.

Definition

DICOM (Digital Imaging and Communications in Medicine) is an open standard for storing and exchanging information in medical image data management.

Function call

Requirement: DICOM must be active.

- Actuate the <Import patient> button in the "Patient database" window. The "Import patient information" window is displayed.

"Import patient information" window

"Select patient" table

Table with personal information which e.g. is already stored in the clinic.

- Actuate a line in the table to select the patient concerned. The table line is color highlighted. The information of the selected patient is displayed in the fields above the table.

Search functions

The fields above the table are input and menu fields which are used for the search functions. An input or menu field is provided for each column of the "Select patient" table:

- "Patient ID"
- "Family name"
- "Given name"
- "Date of birth"
- "Sex"

You can search for several criteria at the same time through simultaneous entries in multiple fields.

You also can search using the wildcard function, i.e. partial inputs and an appended * (asterisk).

<Search> button

When you have entered your search criteria, you can trigger the search function by actuating this button.

If any patient information corresponding to your search criteria is found, it will be displayed in the table.

"Select Patient" table

A table containing personal data which may already be saved at the clinic.

- Activate a line in the table in order to select the corresponding patient. The table line is highlighted in color. The data for the selected patients will be displayed in the fields above the table.

<Import> button

The selected patient files can be imported to your INTRABEAM 600 by actuating the <Import> button.

You still must complete the imported patient information in the INTRABEAM 600 by adding a diagnosis.

<Exit> button

Actuating the <Exit> button returns you to the "Patient Database" window.

Import patient

Search

Show All

44331	Ames	Ames	1942-07-13	Male
42694	Alshen	Alshen-Ingov	1967-05-09	Other
4576	Apone Criga	Apone Criga	1952-01-18	Female
43803	Aaron	Chen	1940-05-09	Female
41468	Ardon	Chen	1911-07-34	Male
42961	Ardia	Arora	1983-09-08	Male
41010	Asper	Bur	2001-10-20	Other
43340	Curea	Uz de	1949-04-14	Other
48230	Fence	Fence	1914-05-04	Female
48825	Egata	Shade	1946-07-25	Other
46290	Fence	Chen	1937-08-16	Other
48840	Emach	Uchi	1986-05-18	Other
41200	Exadib	Exadib	1995-03-03	Other
46990	Spahn	Brake-Cross	1999-11-25	Male
48010	Egata	Shade	1972-02-22	Other
48316	Elata	Olama Sage	1983-07-07	Male
456	Exadib	Exadib	1964-03-17	Other
48241	Exadib	Exadib	1979-11-08	Female
45370	Elu	Yadid	1989-12-21	Male


Page 11 104

5 Menu item <About>

File Tools QA Reports Patient Database About

You can display the basic system information about your INTRABEAM 600, e.g. product name, software and hardware version, serial number, service interval etc. in a window.


- Actuate the menu item <About>;
The "About" window is displayed.

Service interval and service timer

In the "Information" window, i.e. the prescribed service interval and service timer, i.e. the number of days remaining in the current service interval, are displayed (1).

Once the remaining number of days reaches a certain minimum, a warning is displayed during each login to the INTRABEAM software.

- Contact ZEISS Service in time.

Operation of your INTRABEAM 600 will be restricted to service personnel, physicists and the administrator until the service work has been performed. Following an additional password entry, the following users can temporarily continue working with the INTRABEAM 600.

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page. There is no handwriting or other markings on the paper.

Equipment inventory

The "Equipment inventory" is used to record and display the data of all components, i.e. the entire equipment of your INTRABEAM 600. The "Equipment inventory" also is used to monitor the periodically required calibrations and to maintain the equipment data.

Histories of the individual components give you an overview of their individual calibration status and what persons have performed calibrations on them.

Function call

- Actuate the <Equipment inventory> tab;
the "Equipment inventory" window is displayed.

The following groups are provided in the "Equipment inventory":

- 1 "XRS linear accelerator" group
In this group, you administer all (XRS 4) X-ray sources of your INTRABEAM 600; for more information, see Page 142.
- 2 "Equipment" group
In this group, you administer all equipment/components for quality assurance of your INTRABEAM 600; for more information, see Page 148.
- 3 "Applicators" group
In this group, you administer all applicators of your INTRABEAM 600; for more information, see Page 152.

In order to edit the data of a group manually, actuate the corresponding <Edit> button.

Risk due to hardware or data carrier defect

- Always import calibration data prior to a treatment in the context of equipment inventory or by the beginning of quality assurance at the latest.

Version 2.0
G-30-1987-en

Equipment inventory, "XRS linear accelerator" group

In this group, you administer all of the (XRS 4) X-ray sources you have in your INTRABEAM 600.

The data for the XRS 4 that is supplied along with your INTRABEAM 600 are calibration data from ZEISS. These data are installed by ZEISS Service during the initial installation of your INTRABEAM 600.

You can enter additional data manually; see "Importing calibration data" Page 119.

Function call of "XRS inventory"

- Go to the "Equipment inventory" window (see Page 140) and the "XRS linear accelerator" group and actuate the <Edit> button.
The "Edit XRS linear accelerator" window is displayed.

Control functions of the "Edit XRS linear accelerator" window

All of the (XRS 4) X-ray sources registered in your INTRABEAM 600 are listed in a table (1) in this window.

Individual data (columns of the table) can be edited. Corresponding buttons and input fields are provided above the list for this purpose.

1 Table, selecting the XRS 4 (table line)

Actuate the corresponding line directly in the table in order to edit the data of a specific CRS 4.

The selected line is color highlighted. The editable data in the selected line are displayed above the table in the corresponding menu buttons and input fields, where they can be edited.


The calibration data of the XRS 4 concerned are displayed in a separate window. To do this, actuate the table line concerned with a double-click (of the mouse).

2 "XRS name" input field.

displays the name of the currently selected XRS 4. This name was internally assigned for this XRS 4 in your INTRABEAM 600. You can edit this name.

3 <Status> menu button

The status of the currently selected XRS 4 is displayed in the button. When the button is actuated, a menu with the menu items <Active> and <Inactive> opens.

4 <Validity> menu button

The calibration validity period (of the selected XRS 4) is displayed in the button. When the button is actuated a menu opens. You can select a different validity period.

5 "Notes" input field

The currently saved notes on the selected XRS 4 are displayed. These notes can be edited in this field.

6 <Apply> button

The current settings and inputs for the menu buttons and input fields are accepted for the selected XRS 4. After actuating the <Accept> button, you can change to a different XRS 4 (table line).

7 <Delete> button

The XRS 4 highlighted in the table (line highlighted in the table) is deleted from the table and is no longer available.

8 <Exit> button

Rejects all current changes which have not yet been saved by actuating the <Accept> button and returns the program to the parent window, i.e. "Equipment inventory".

XRS name	XRS S/N	Voltage [kV]	Calibrated by	Cal date	Validity	Verified by	Status	DR [Gy/min]	Check date	Notes
2nd	507528	40.0	226 / 226	22-Jan-16	Annual	Service	Inactive	2.844	22-Jan-16	
	507542	40.0	226 / 226	21-Jan-16	Annual	Service	Inactive	2.937	21-Jan-16	
1st	507542	50.0	226 / 226	21-Jan-16	Annual	Service	Active	3.805	21-Jan-16	

Columns of the "XRS inventory" table

9 XRS name

Name of the XRS 4, individual name in your INTRABEAM 600

10 XRS S/N

Serial number of the XRS 4, specified in the calibration file and on the housing of the XRS 4 concerned.

11 Voltage kV

50 kV, 40 kV: The voltage determines the maximum power of the X-ray radiation

12 Calibrated by

Name entered automatically

13 Cal. date

entered automatically

14 Validity

You must select the period of validity for each XRS 4 manually with the "Validity" button (4).

- Observe the applicable legal regulations for the calibration intervals of the XRS 4; we recommend annual calibration.
- Observe the calibration regulations on Page 254.

When the calibration interval expires, the XRS 4 concerned is automatically set to "Inactive".

The calibration interval of the XRS 4 starts with the initial startup of your INTRABEAM 600. Since several months may possibly pass between the calibration and the initial startup, "Biannual" is entered in this field on initial delivery. Nevertheless, ZEISS will automatically call your attention to the expiry of the calibration interval within one year following initial startup.

- Reset the calibration interval to "Annual" when the first calibration interval expires.

15 Verified by

Name of the person who imported and accepted the calibration file.

16 Status

The status of the respective XRS 4 is displayed in the "Status" column.

"Active": The respective XRS 4 can be used for verification.

"Inactive": The respective XRS 4 is not available for verification.

Using the <Status> menu button (3), you can switch the respective XRS 4 from "Inactive" to "Active".

If you set an XRS 4 with an expired calibration to "Active", you must confirm this with your password, i.e. all further operation will be performed on your own personal responsibility.

You must repeat this activation following each software restart.

17 DR [Gy/min]

Dose rate, result of quality assurance "Output check".

18 Check date

Date of last quality assurance.

19 Notes

Individual comments can be stored for each XRS 4 using the input field (5).

XRS linear accelerator																			
XRS number		Voltage		Calibration by		Cal date		Validity		Verified by		Status		Left [kGy/min]		Check date		Notes	
		Inactive		226		Annual						Apply				Delete			
XRS number	XRS S/N	Voltage [kV]	Calibration by	Cal date	Validity	Verified by	Status	Left [kGy/min]	Check date	Notes									
2nd	507538	50.0	226 / 226	22-Jan-16	Annual	Service	Active	3.589	22-Jan-16										
	507542	40.0	226 / 226	21-Jan-16	Annual	Service	Inactive	2.917	21-Jan-16										
1st	507547	50.0	226 / 226	21-Jan-16	Annual	Service	Active	3.005	21-Jan-16										
9	10	11	12	13	14	15	16	17	18	19									

Equipment summary

Quality assurance

Treatment prescription

Treatment delivery

Treatment monitoring

Treatment review

Calibration

Recalibrating the XRS 4

We recommend having your XRS 4 calibrated annually. Your XRS 4 is calibrated by ZEISS at the manufacturer's plant.

During calibration, the characteristic curve, dose rate [Gy/min] in relation to the penetration depth [mm] of the respective XRS 4 is recorded at 40 kV and at 50 kV. This takes place in a water phantom, i.e. in a special container filled with water.

Your recalibrated XRS 4 will be returned to you by ZEISS with a current calibration file. You must copy this calibration file to the software in your computer, see section on "Importing calibration data" Page 119.

- The old calibration file is automatically moved to the archive folder.

Saving the calibration file

The software detects the availability of new calibration data. The software displays the following prompt during the next start: "Please confirm the XRS calibration data".

If you confirm this prompt with <Accept>, the new calibration data will be used immediately.

The same window also will be displayed later on if you select an XRS 4 (line in the table) with a double-click of the mouse. In this case, the window is displayed only for inspection.

Quality assurance

Treatment prescription

Treatment results

Calibration

01 Calibration verification

RSC

1.01.1.1.01.01

50 mV

21-Jan-16


VNA 10

40 uA


226 / 226

907542

TARGET

Calibration: 10

Accept

Page 10

Equipment inventory, "Equipment" group

In this Equipment inventory group, you administer equipment especially for the quality assurance of your INTRABEAM 600:

- CC 600 control unit
- IC ionization chamber
- PAICH (probe adjuster and ionization chamber holder)
- NC32 Floor Stand
- UNIDOS E dosimeter
- PDA (photo diode array)

The data for the components supplied along with your INTRABEAM 600 are entered by ZEISS Service when first setting up your INTRABEAM 600 in this group. You can enter additional data manually.

"Equipment" function call

- In the "Equipment inventory" window (see Page 140), actuate the <Edit> button belonging to the "Equipment" group.
The "Edit equipment" window is displayed.

Control functions of the "Edit equipment" window

All of the equipment/components registered in your INTRABEAM 600 are listed in a table (1) in this window.

All columns of this table concerning further components of the INTRABEAM 600 can be edited here. The appropriate buttons, menu buttons and input fields are provided for this purpose above the table.

1 Table, selecting the component (table line)

Tap directly on the corresponding row in the table to edit the data for a particular component.

The selected line is color highlighted. The editable data in this line are displayed in the corresponding menu buttons and input fields above the table, where you can edit them.

2 <Apply> button

The current settings and inputs for the menu buttons and input fields are accepted for the selected component. After actuating the <Accept> button, you can change to a different component.

**3** <New> button

Enter a new component in the table, i.e. a new line is created in the table.

4 <Delete> button

Delete currently selected component, i.e. delete the line currently selected in the table.

5 <Type> menu button

The currently selected component is displayed. To make changes or to make a new entry, tap the menu button to open the relevant menu and select a menu item.

6 <Model> menu button

The model of the currently selected component is displayed. To make changes or to make a new entry, tap the menu button to open the relevant menu and select a model variant.

The screenshot shows the application interface with various fields and buttons. Numbered callouts (1-16) point to specific elements:

- 1: Apply button
- 2: New button
- 3: Delete button
- 4: Calibration factor field
- 5: Equipment Type dropdown menu
- 6: Model dropdown menu
- 7: S/N field
- 8: Calibration factor field
- 9: Calibration unit dropdown menu
- 10: Verified by field
- 11: Calibrated by field
- 12: Cal date (DD.MM.YY) field
- 13: Status dropdown menu
- 14: Notes field
- 15: Notes field (continued)
- 16: Exit button

Type	Model	S/N	Cal factor	Cal unit	Cal date	Validity	Cal by	Verified by	Status	Notes
Control Console	CC 600	1	0.000		14-Apr-16	Permanent			Active	
PDA	PDA	507315	0.000		14-Apr-16	Permanent	CZ	CZ	Active	
PAICH	PAICH	507401	-27.465	°C	14-Apr-16	Permanent	CZ	CZ	Active	
Ion Chamber	TN23342A	002013	1.081F+9	Gy/C	29-Oct-15	Annual	PTW	CZ	Active	
Dosimeter	UNIDOS-E	001025	1.000		14-Apr-16	Annual			Active	
Floor Stand	NEW	000001	0.000		25-Apr-16	Biennial			Inactive	

7 <S/N> input field

The INTRABEAM 600 requires the serial number of the respective component for communication and logging. The serial number of the component can be entered or edited here.

8 <Calibration factor> input field

The calibration company supplies the calibration factor to you along with the calibration certificate of the component concerned.

For example, the calibration factor for the IC ionization chamber: The calibration factor with the unit of measure Gy/C is used to measure the dose rate during the "Output check", see Page 182.

PAICH data: The data of the PAICH probe adjuster and ionization chamber holder come from the measured result of the calibration procedure "Calibrating the temperature sensor", <Thermometer>, see Page 267. These data are automatically entered following the calibration procedure.

9 <Calibration unit> input field

Unit of measure of calibration factor, for example Gy/C. The calibration company supplies the calibration unit to you along with the calibration certificate of the component concerned.

10 <Verified by> input field

Name of checker who checked the newly calibrated component and accepted it for the INTRABEAM 600.

11 <Calibrated by> input field

The calibration company supplies the name of the calibrator to you along with the calibration certificate of the component concerned.

12 <Cal. date (DD MM YY)> input field

The calibration company supplies the date of calibration to you along with the calibration certificate of the component concerned.

13 <Status> menu button

The status of the currently selected component is displayed in the button. When the button is actuated, a menu with the menu items <Active> and <Inactive> opens.

"Active" The component is available for further procedures.

"Inactive" The component is no longer available.

Only one component type can be active at a time!

**14 <Validity> menu button**

Observe the applicable regulations for the different instruments.

15 <Notes> input field

Individual comments can be stored and edited using the input field.

16 <Exit> button

Rejects all current changes which have not yet been saved by actuating the <Accept> button and returns the program to the parent window, i.e. "Equipment inventory".

Equipment inventory

Equipment Type: Control Console, Model: CC 600, S/N: 1, Calibration factor: 0.0000E+0, Calibration date: , Validity: Permanent, Status: Active, Cal date (DD.MM.YY): 14.04.16, Calibrated by: , Verified by: , Notes:

Buttons: Apply, New, Delete

Type	Model	S/N	Cal factor	Cal unit	Cal date	Validity	Cal by	Verified by	Status	Notes
Control Console	CC 600	1	0.000		14-Apr-16	Permanent			Active	
PDA	PDA	507315	0.000		14-Apr-16	Permanent	CZ	CZ	Active	
PAICH	PAICH	507401	-27.465	°C	14-Apr-16	Permanent	CZ	CZ	Active	
Ion Chamber	TN23342A	002013	1.081E+9	Gy/C	29-Okt-15	Annual	PTW	CZ	Active	
Dosimeter	UNIDOS-E	001025	1.000		14-Apr-16	Annual			Active	
Floor Stand	NEW	000001	0.000		25-Apr-16	Biennial			Inactive	

Buttons: Exit

Equipment inventory, "Applicators" group

In this group of Equipment inventory, you administer all of the applicators you have in your INTRABEAM 600.

The data for the applicators that are supplied along with your INTRABEAM 600 are calibration data from ZEISS. These data are installed by ZEISS Service during the initial installation of your INTRABEAM 600.

You can enter additional data manually; see "Importing calibration data" Page 119.

"Applicators" function call

- In the "Equipment inventory" window (see Page 140), actuate the <Edit> button belonging to the "Equipment" group.
The "Edit applicators" window is displayed.

All of the applicators registered in your INTRABEAM 600 are listed in a table (1) in this window.

Individual data (columns of the table) can be edited. Corresponding buttons and input fields are provided above the list for this purpose.

1 Table, selecting the applicator (table line)

Tap directly on the corresponding row in the table to edit the data for a particular applicator.

The selected line is color highlighted. The editable data in this line are displayed in the corresponding menu buttons and input fields above the table, where you can edit them.

The calibration data of the applicator concerned are displayed in a separate window. To do this, actuate the table line concerned with a double-click (of the mouse).

2 <Status> menu button

The status of the currently selected applicator is displayed in the button. When the button is actuated, a menu with the menu items <Active> and <Inactive> opens.

"Active": The applicator concerned can be used for a further treatment.

"Inactive": The applicator concerned is blocked for a further treatment.

If you set an applicator which has reached the maximum permissible number of sterilizations to "Active", you must confirm this with your password, i.e. all further operation will be performed on your own personal responsibility.

Electronic dosimetry

Quality assurance

Treatment prescription

Treatment plan

Treatment advice

Calibration

**WARNING****Permissible number of processings exceeded**

Once an applicator's permissible number of processings has been exceeded, its cleanliness, sterility and biocompatibility no longer can be guaranteed. The functioning of the applicator no longer can be ensured.

Defective applicators may result in the emission of a radiation dose which does not correspond to the planned dose. The consequences are underexposure or radiation damage due to overexposure. If the patient or user comes into contact with an unsterile or no longer biocompatible applicator, the consequence of this is allergic reactions and infections.

- Ensure that all processing cycles are counted and that no applicator which has undergone more than the permissible number of processing cycles is used. For the maximum permissible number of processing cycles, see the table below.

Applicator type	Maximum permissible number of processing cycles
Spherical applicator	100
Surface applicator	50
Flat applicator	50

1 2 3 4 5 6 7 8

Applicators

Status: Inactive

Usage: 7

Notes: n/a

Apply

Delete

Show All

Show All

Type	Size (cm)	SyRa	Calibrated by	Calibration date	Verified by	Usage	Status	Notes
Spherical	1.5	A61065	PL	22-Okt-15	Service	7/100	Inactive	n/a
Spherical	1.5	A63005	zb	08-Okt-15	Service	0/100	Inactive	n/a
Spherical	2.0	B63001	zb	08-Okt-15	Service	1/100	Inactive	n/a
Spherical	2.0	B63002	zb	08-Okt-15	Service	0/100	Inactive	n/a
Spherical	5.0	H63002	zb	08-Okt-15	Service	1/100	Inactive	n/a
Spherical	5.0	H63001	zb	08-Okt-15	Service	0/100	Inactive	n/a
Needle	0.5	N20010	zb	08-Okt-15	Physicist	0/1	Inactive	n/a
Needle	0.5	N20011	zb	08-Okt-15	Physicist	0/1	Inactive	n/a
Needle	0.5	N11889	PL	21-Okt-15	Service	0/1	Inactive	n/a
Balloon	3.5	310002	PL	29-Jan-14	Service	0/0	Inactive	n/a

Exit 9

Equipment inventory	Quality assurance	Treatment prescription	Treatment delivery	Treatment review	Calibration
---------------------	-------------------	------------------------	--------------------	------------------	-------------

3 <Usage> input field

The number of previous usages (sterilization cycles) is displayed. An applicator is automatically incremented for each usage. When the prescribed maximum number of usages has been reached, the applicator is automatically set to "Inactive".

You may sterilize an applicator only a limited number of times. Then you are no longer allowed to use it.

For radiotherapy treatment, several applicators of varying size are often provided in the operating room. All these applicators are then regarded as "contaminated" and must be reprocessed.

Only an applicator selected in the context of "treatment prescription" is automatically incremented by the INTRABEAM 600 as soon as the treatment begins. You must document the processing of all other applicators manually. The <Usage> input field is used for this purpose.

4 <Notes> input field

The currently saved notes on the selected applicator are displayed. These notes can be edited in this field.

5 <Apply> button

The current settings and inputs for the menu buttons and input fields are accepted for the selected applicator. After actuating the <Accept> button, you can change to a different component.

6 <Delete> button

The applicator highlighted in the table (line highlighted in the table) is deleted from the table.

7 <Show All> menu button, "Status" display filter

When the button is actuated, a menu with the filter criteria for the display opens in the table. With this menu button, you can make selections under the menu items <Active> and <Inactive>. Once you have selected a menu item, only the applicators that meet the filter criterion will be displayed in the table.

8 <Show All> menu button, "Applicator Type" display filter

When the button is actuated, a menu with the filter criteria for the display opens in the table. You can use this menu button to select one of the applicator types. Once you have selected a menu item, only the applicators that meet the filter criterion will be displayed in the table.

9 <Exit> button

Rejects all current changes which have not yet been saved by actuating the <Accept> button and returns the program to the parent window, i.e. "Equipment inventory".

Type	Size [cm]	S/N	Calibrated by	Calibration date	Verified by	Usage	Status	Notes
Spherical	1.5	A61065	PL	22-Oct-15	Service	7/100	Inactive	n/a
Spherical	1.5	A63005	zb	08-Oct-15	Service	0/100	Inactive	n/a
Spherical	2.0	B63001	zb	08-Oct-15	Service	1/100	Inactive	n/a
Spherical	2.0	B63002	zb	08-Oct-15	Service	0/100	Inactive	n/a
Spherical	5.0	H63002	zb	08-Oct-15	Service	1/100	Inactive	n/a
Spherical	5.0	H63001	zb	08-Oct-15	Service	0/100	Inactive	n/a
Needle	0.5	N20010	zb	08-Oct-15	Physicist	0/1	Inactive	n/a
Needle	0.5	N20011	zb	08-Oct-15	Physicist	0/1	Inactive	n/a
Needle	0.5	N11889	PL	21-Oct-15	Service	0/1	Inactive	n/a
Balloon	3.5	310002	PL	29-Jan-14	Service	0/0	Inactive	n/a

Buttons: Apply, Delete, Show All, Show All, Exit

Columns of the "Applicator inventory" table

10 Type

Type of the respective applicator, from calibration data, not editable

11 Size [cm]

Size of the respective applicator, from calibration data, not editable

12 S/N

Serial number of the respective applicator, from calibration data, not editable. The INTRABEAM 600 requires the serial number of the respective applicator for communication and logging.

13 Calibrated by

Name specification, from calibration data, not editable

14 Calibration date

Date specification, from calibration data, not editable

15 Verified by

Name specification, from calibration data, not editable

16 Usage

The number of usages and the number of sterilization cycles therefore required is automatically incremented and displayed here for those applicators which are actually used.

In the case of applicators which become unsterile without actually being used, you must increment the number of usages manually.

When the prescribed maximum number of usages has been reached, the applicator is automatically set to "Inactive".

You may sterilize an applicator only a limited number of times. Then you are no longer allowed to use it.

17 Status

Current status according to the number of usages, editable.

18 Notes

Any currently saved notes concerning the applicator in question will be displayed. These notes can be edited in the "Notes" input field.

10	11	12	13	14	15	16	17	18
Item	Size	Part	Material	Lot	Service	Exp.	Status	Notes
Spherical	1.5	A63005	st	08-08-15	Service	0/100	Inactive	n/a
Spherical	2.0	B63001	st	08-08-15	Service	1/100	Inactive	n/a
Spherical	2.0	B63002	st	08-08-15	Service	0/100	Inactive	n/a
Spherical	3.0	H63002	st	08-08-15	Service	1/100	Inactive	n/a
Spherical	3.0	H63001	st	08-08-15	Service	0/100	Inactive	n/a
Needle	0.5	N20010	st	08-08-15	Physicist	0/1	Inactive	n/a
Needle	0.5	N20011	st	08-08-15	Physicist	0/1	Inactive	n/a
Needle	0.5	N11889	Pl	21-08-15	Service	0/1	Inactive	n/a
Balloon	3.5	310002	Pl	20 Jan 14	Service	0/0	Inactive	n/a

Equipment inventory

Quality assurance

Treatment prescription

Treatment plans

Calculations

Quality assurance

Prior to each treatment, you must perform the prescribed verifications on the XRS 4 to be used with the help of "Quality assurance". Certain verification processes are required for quality assurance which you must perform prior to each treatment. Other verification processes are recommended.

For the INTRABEAM 600, verification means the metrological confirmation and verification that the X-ray source (XRS 4) meets the quality requirements for radiation emission.

For quality assurance, you must at least perform the "Isotropy check" and the "Output check" successfully.

The following procedures are used for quality assurance:

- Probe adjuster (XRS tip adjustment)
The straightness and centering of the X-ray probe can be checked and corrected with the "PAICH" probe adjuster and ionization chamber holder. You must perform the XRS tip adjustment if the isotropy check was not successfully completed or if the X-ray probe could possibly be bent.
- Dynamic offset
The alignment of the electron beam to the mechanical center of the X-ray probe can be adjusted via beam deflection using the deflection coils. For this purpose, the direction of the electron beam can be checked with the "PDA" photo diode array.
A mechanical alignment of the x-ray probe (XRS tip alignment) must always be followed by an electronic alignment (dynamic offset).
- Isotropy check
You must always check the geometry of the spherical radiation emission on the probe tip. You require the PDA photo diode array for the isotropy measurement.


Procedure	
Probe adjuster	Recommended
Dynamic offsets	Recommended
Isotropy check	Mandatory
Output check	Mandatory

– **Output check**

You measure the constancy of the dose rate of the XRS 4 using the IC ionization chamber and the dosimeter.

The dose rate measured by you is compared with a specific expected value from the calibration data. A correction factor between the measured dose rate and the data in the calibration file is calculated and taken into account later on in the following treatment planning. The prerequisite is a successfully completed isotropy check.

Quality assurance workflow

Equipment inventory

Quality assurance

Treatment prescriptions

Treatment review

Calibration

Performing quality assurance

NOTE

Danger in case of nickel allergy

If a user with a nickel allergy has skin contact with the housing of an XRS 4 or an inspection instrument (PAICH, PDA), allergic reactions may result.

- Always wear suitable gloves and protective clothing if necessary.

Function call

- Actuate the <Quality Assurance> tab.
The "Quality Assurance" window is displayed.

NOTE

If the interlock surfaces are contaminated, it is possible that the interlock will not be closed and the INTRABEAM 600 will not generate any X-ray radiation.

- Ensure (visually) that the interlock surfaces are clean prior to every use.

1 Select XRS 4

- Actuate the <S/N> keypad.
A menu listing all XRS 4s which are set to "Active" in the equipment inventory of your INTRABEAM 600 (Page 142) opens. An XRS 4 is specified by its serial number and voltage (kV) for each menu item.
- Select the XRS 4 that you would like to use for treatment later on.

2 Validity of verification

If the "QA completed" indicator lamp lights up green, this means the verification was successful. The verification of the XRS 4 is now valid for a maximum of 36 hours. Starting from now, a countdown of 36 hours is performed. The remaining time for which the verification will remain valid is displayed under "Active for" next to the indicator lamp. The time display is updated every time an XRS 4 and the associated voltage are selected.

The verification loses its validity after 36 hours or due to completion of the treatment and "0.0 h" is displayed next to the indicator lamp under "Active for".

If the isotropy test is successfully completed again within 36 hours following successful verification, this does not automatically cause the constancy test to be set to "Required". The "Constancy test" indicator light remains green.

3 "Mandatory" or "Recommended" quality assurance procedures

Procedure


Probe adjuster	Recommended
Dynamic offsets	Recommended
Isotropy check	Mandatory
Output check	Mandatory

"Mandatory"

For procedures marked "Mandatory", verification of the measured values concerned is mandatory, i.e. a treatment can be performed only if these requirements are met (Color code "green").

"Recommended"


For procedures marked "Recommended", verification of the measured values concerned is required only if the isotropy check for the XRS 4 was not successful. Otherwise these procedures are optional.

Color code for quality assurance procedures

For the selected XRS 4, the status of each individual procedure is displayed by the colored indicator light in the respective button.

Color yellow The selected XRS 4 is ready in the system for this verification procedure

Color green The selected XRS 4 has been successfully verified for procedure.


Color red Error; the selected XRS 4 cannot be used until it has passed this verification, see "Workflow" Page 104.

4 "QA completed" display (4)

When you have successfully completed the quality assurance of the respective XRS 4, the "QA completed" indicator light lights up green. If the quality assurance steps were not performed correctly, the INTRABEAM 600 will not allow irradiation with the respective XRS 4.

QA Reports

The QA Reports of your INTRABEAM 600 can be opened via the menu bar in the <QA Reports> menu, see Page 133:

<QA Reports> / <Open>
<QA Reports> / <Open Latest>

The QA Reports of your INTRABEAM 600 can be exported using the "File Browser" window, see Page 120:

- Actuate the <Data Management> button in the INTRABEAM Start Manager; the "File Browser" window is displayed.
- Select the folder "QA Reports" in the "File Browser" window.

Equipment inventory Quality assurance Treatment prescription Treatment review Calibration

1 2 3 4

Serial number: 397542 50.0 kV

Serial photoresistor ID: 0.0 hr

Serial photoresistor ID: 0.0 hr

Procedure:

- * Recommended
- * Recommended
- * Mandatory
- * Mandatory

Radial beam position:

0.00 mm 0.00 mm

0.00 mm 0.00 mm

Press Start to run PDA Source Check at 50.0 kV, 40 uA or Stop to terminate

Start

Stop

Equipment inventory

Quality assurance

Treatment prescription

Treatment device

Calibration

Proper handling of components during quality assurance


WARNING

Risk of damage when inserting the X-ray probe in an inspection component

The X-ray probe may be mechanically damaged or rendered useless by careless insertion in one of the inspection components (PAICH or PDA). This could possibly lead to false irradiation.

- Always use the V-guide to insert the X-ray probe in one of the inspection components.



WARNING

Danger if V-guide is used with incorrect orientation

If, while inserting the X-ray probe in an inspection component (PAICH or PDA), the V-guide is used with an incorrect orientation, the X-ray probe may be mechanically damaged and thus rendered useless. This can possibly lead to false irradiation.

- Note that one side of the V-guide is designed especially so that it can receive the XRS 4, while the other side can receive the verification instruments. The V-guide is marked with a corresponding engraving.

- To set up the XRS 4 upright with the PDA or PAICH, place the XRS 4 into the XRS mount with an exactly vertical orientation; this prevents it from falling over.
- Connect the PDA / PAICH and the XRS 4 to their respective cables.

Equipment inventory

Quality assurance

Treatment prescription

Treatment review

Calibration

NOTE**Crushing hazard when inserting the X-ray probe in an inspection component**

If you carelessly insert the X-ray probe in one of the inspection components (PAICH or PDA), your fingers may get caught between the two housings.

- Always use the V-guide to insert the X-ray probe in one of the inspection components.

**CAUTION****Danger due to exposed beryllium (Be)**

The probe tip of the XRS 4 contains beryllium (Be), and therefore has a protective coating.

The protective coating of the probe tip may, for example, be damaged if the probe tip is bumped while being inserted in the inspection component. If the user then touches the damaged probe tip, he may contaminate his skin and clothing with beryllium.

- Always avoid direct skin contact with X-ray probes.
First aid measure: Following skin contact with beryllium, wash with plenty of water. Remove contaminated clothing.
- If a probe tip is damaged, return the XRS 4 to ZEISS immediately.

NOTE**Risk of damage when removing the X-ray probe from the inspection component**

If the X-ray probe is removed from one of the inspection components (PAICH or PDA) carelessly and without a guide, this could cause damage to the X-ray probe.

- Always use the V-guide to remove the X-ray probe from the inspection component.

Video instruction, video window

When you start a quality assurance or calibration procedure, the steps to be performed are displayed in video instructions with text overlay.

The software internally checks whether all of the required inspection components are connected and in position. If this is not the case, the selected quality assurance or calibration procedure cannot be continued.

Notes on the video window

The example shown here is the video window on the quality assurance procedure "Checking the XRS 4 with the PDA".

If the latest quality assurance or calibration procedure is used, the corresponding video sequence will be displayed in field (1).

Descriptive text on the current quality assurance or calibration procedure is displayed in field (2).

The inspection component which detected the optical interlock (in this example, the PDA) is displayed in field (3).

The status displays (4) indicate the interlock statuses of the cable connections. If the corresponding status display (4) lights up green, the required conditions have been fulfilled, here for example:

- XRS cable connected
 - PDA cable connected
 - Status of external interlock
 - With the <Cancel> button (5), you can cancel the current quality assurance or calibration procedure.
 - Actuate the <Continue> button (5) to continue the current quality assurance or calibration procedure.
- When the status display is not lit or the interlock status hasn't been fulfilled, the <Continue> button (5) is not active.

"Video information" window

1 Video display field

The display field provides you with instructions for use regarding the currently selected screen mask in the form of a video sequence.

2 Text display field


This display field provides you with instructions for use regarding the currently selected screen mask in text form.

3 Specification of current component,
i.e. status display of optical interlock (here in the example "PAICH").
The picture below shows the video window concerning the quality assurance procedure "Checking the XRS 4 with the PDA".

4 Status displays
If the corresponding status display lights up green, the required conditions have been fulfilled.

5 <Cancel> button
You can use this button to cancel the selected quality assurance procedure.

6 <Continue> button
You can continue the procedure only if the status displays are lit and the interlock condition has been fulfilled.

Equipment inventory

Quality assurance

Treatment prescription

Treatment review

Calibration

XRS tip adjustment

The "XRS tip adjustment" is used to check and correct the straightness and centering of the X-ray probe with the probe adjuster and ionization chamber holder (PAICH).

You absolutely must perform the "XRS tip adjustment" if the "Isotropy check" (Page 176) was not successfully completed or if it is suspected that the X-ray probe could possibly be bent. If the "Isotropy check" was successful, the "XRS tip adjustment" is optional.



If the XRS 4 already was successfully checked and you have subsequently centered the X-ray probe with the "XRS tip adjustment", it is mandatory to repeat the "Dynamic offset", "Isotropy check" and "Output check".

- If you have a PAICH with a micrometer screw, you must ensure that the micrometer screw cannot damage the X-ray probe: i.e. before inserting the X-ray probe in the PAICH, you must unscrew the micrometer screw on the PAICH until the marking "3" becomes visible on the micrometer scale.

The proper procedure for inserting the X-ray probe in the PAICH using the V-guide and then setting up the XRS 4 in the XRS mount with the PAICH mounted vertically is described in the section "Proper handling of components during quality assurance", see Page 164.

- Insert the X-ray probe in the PAICH using the V-guide.

- To set up the XRS 4 upright with the PAICH, place the XRS 4 into the XRS mount with an exactly vertical orientation; this prevents it from falling over.
- Connect the PAICH and the XRS 4 with their respective cables.

Equipment inventory

Quality assurance

Treatment prescription

Treatment review

Calibration

Acclimitization

Place the PAICH in any starting position. Let the setup acclimatize for at least one minute before starting the measurements.

Procedure

Probe adjuster	Recommended
Dynamic offsets	Recommended
Isotropy check	Mandatory
Output check	Mandatory

Function call

- Actuate the <Probe adjuster> button in the "Quality assurance" window.

Start of video instructions

As soon as one of the quality assurance or calibration procedures has been started, i.e. one of the corresponding buttons has been actuated, the related video instructions automatically start.

The steps to be performed are displayed in the respective video instructions as animated graphics with text, see Page 166.

The software internally checks whether all of the required inspection instruments are connected and in position. If this is not the case, the procedure just started cannot be continued.

Operational sequence for XRS tip adjustment

- Mount the PAICH cable in the holder on the monitor.
- Follow the instructions provided (6) in the "XRS tip adjustment" window.
- Actuate the <Zero> button (4).
The displays of the deviations are reset to zero.
- Slowly rotate the PAICH slightly more than 360°:
 - For this purpose, grasp the PAICH by its base.
 - Do not change your grip too often when rotating the PAICH.

If you do not rotate the PAICH by at least 360° during a measurement, that measurement will be incomplete.

If a deviation of less than 0.1 mm is displayed, the straightness of the X-ray probe is within tolerance.

- You can terminate the XRS tip adjustment by actuating the <Done> button.

Under "Total deviation" the sum of the measured deviations is displayed.

If a deviation of more than 0.1 is displayed, the X-ray probe requires adjustment:

- Place the PAICH in any starting position. Let the setup adapt to room temperature for at least one minute.
- Rotate the PAICH slowly through 360° so that the complete deviation of the XRS tip can be measured; i.e. turn the PAICH until the deviation bars display their highest positive value.
- Now calculate the target value to be set with the pusher released: Current deviation (left-hand value) minus one half of total deviation.
- Carefully press the pusher until the calculated value appears on the display of the current deviation (2) when the pusher is released.

Experience has shown that the pusher has to be pressed and released several times while observing the bar diagrams until the target value is reached.

- Let go of the pusher and rotate the PAICH back to its starting position.
- Actuate the <Zero> button (4) to reset the deviation displays (3) to 0 mm.
- Repeat the previous steps until the total deviation is less than 0.10 mm, recommendation: approx. 0.07 mm.
- Actuate the <Done> button (5) to accept the result of the XRS tip adjustment and exit the routine.

NOTE

Failure due to excessive eccentricity of probe tip

The eccentricity of the probe tip must be less than 0.10 mm as a prerequisite for enabling the performance of further inspection steps on the XRS 4 concerned. If the eccentricity remains above 0.10 mm, you must repeat the XRS tip adjustment, since this inspection step otherwise has not been passed. You can carry out further steps ("Dynamic offset", Page 173) of the quality assurance process only after successfully completing the XRS tip adjustment.


- Call ZEISS Service if it is not possible to center the probe tip at less than 0.10 mm.

The proper procedure for removing the x-ray probe from the verification components (PAICH / PDA) is described in the section "Proper handling of components during quality assurance", see Page 165.

Legend for "XRS tip adjustment" window

- 1 <Probe adjuster> button
- 2 Deviation bars
- 3 Display of deviations
- 4 <Zero> button
- 5 <Done> button
- 6 "Instruction" display field

[Faint, illegible handwriting visible through the paper]

Equipment inventory

Quality assurance

Treatment prescription

Treatment review

Calibration

Dynamic Offsets (adjusting the beam deflection)

When adjusting the beam deflection, you set the direction of the electron beam to the mechanical center of the X-ray probe with the help of the deflection coils. If the beam deflection is not adjusted, the electron beam may no longer strike the target properly following a realignment of the X-ray probe (tip adjustment) (For more information on the "Target", see Page 74). In this case, the verification steps no longer can be successfully completed.



A mechanical alignment of the X-ray probe, i.e. centering of the X-ray probe, (function: XRS tip adjustment) always must be followed by an electronic alignment of the X-ray probe (function: dynamic offset).

- Always check the beam deflection after centering the X-ray probe (XRS tip adjustment) or importing a new calibration file from a recalibration. Otherwise, this inspection is optional.

The proper procedure for inserting the X-ray probe in the PDA using the V-guide and then setting up the XRS 4 vertically in the X-guide with the PDA mounted is described in the section "Proper handling of components during quality assurance", see Page 164.

- Insert the X-ray probe in the PDA (photo diode array) with the help of the V-guide.

- To set up the XRS 4 upright with the PDA, place the XRS 4 into the XRS mount with an exactly vertical orientation; this prevents it from falling over.
- Connect the PDA and the XRS 4 with their respective cables.
- Align the +X and +Y markings on the XRS 4 and the PDA so that they coincide.

Procedure

Probe adjuster	Recommended
Dynamic offsets	Recommended
Isotropy check	Mandatory
Output check	Mandatory

Function call

- Actuate the <Dynamic offsets> button in the "Quality assurance" window.

Start of video instructions

As soon as one of the quality assurance or calibration procedures has been started, i.e. one of the corresponding buttons has been actuated, the related video instructions automatically start.

The steps to be performed are displayed in the respective video instructions as animated graphics with text, see Page 166.

The software internally checks whether all of the required inspection instruments are connected and in position. If this is not the case, the procedure just started cannot be continued.

- Press <START>. The adjustment is automatically performed for the selected XRS X-ray source at a certain high voltage (XRS@kV). This process takes a few minutes.

You must perform this procedure separately for each additional acceleration voltage.

If the beam deflection is successfully adjusted, the color code display turns "green" in the "Dynamic offset" window, i.e. then you can carry out further inspection steps (isotropy check, Page 176).

If the beam deflection cannot be successfully adjusted, an error message will appear to indicate the reason why. In this case, please contact ZEISS Service. If the color code display for PDA and PAICH already has turned "green", you must repeat these.


The proper procedure for removing the x-ray probe from the verification components (PAICH / PDA) is described in the section "Correct handling of components during quality assurance", see Page 165.

Legend for "Dynamic offset" window

- 1 <Dynamic offsets> button for function call
- 2 Graphic representation of beam deflection from the center
- 3 Display of "X adjustment" and "Y adjustment" offset values
- 4 Text display of current inspection steps
- 5 <START> button
- 6 <STOP> button

Isotropy check

You must always perform this check!

The geometry of the spherical radiation emission on the probe tip can be checked with the "Isotropy check". The isotropy check and isotropy adjustment are performed at maximum beam current and the acceleration voltage selected for the respective XRS 4.

The INTRABEAM 600 permits treatment only using a setting of the acceleration voltage for which the "isotropy check" has been successfully completed.

Performing an "isotropy check" invalidates the result of the "output check" (XRS 4 with PAICH), which is a lower-ranking quality assurance test; accordingly, this test must be repeated before treatment can begin.

You require a PDA (photo diode array) for the isotropy adjustment. Operation and display occur in the "Isotropy check" window.



- Perform the "Isotropy check" procedure on an XRS 4 with all of the acceleration voltages in question several times before going on to the remaining inspection tests. This prevents having to repeat these tests in case an acceleration voltage other than the one you used while performing the isotropy adjustment is required for a treatment.

The proper procedure for inserting the X-ray probe in the PDA using the V-guide and then setting up the XRS 4 vertically in the XRS mount with the PDA mounted is described in the section "Proper handling of components during quality assurance", see Page 164.

- Insert the X-ray probe in the PDA with the help of the V-guide.

- To set up the XRS 4 upright with the PDA, place the XRS 4 into the XRS mount with an exactly vertical orientation; this prevents it from falling over.
- Connect the PDA and the XRS 4 with their respective cables.

Isotropy adjustment (adjustment of the radiation field on the probe tip)

The adjustment of the radiation field on the probe tip essentially comprises two steps:

- Setting the operating parameters (acceleration voltage, beam current and running time)
- Automatic determination of the isotropy of the radiation field on the probe tip.

Procedure

☒ Probe adjuster	Recommended
☒ Dynamic offsets	Recommended
☒ Isotropy check	Mandatory
☒ Output check	Mandatory

Function call

- Actuate the <Isotropy check> button in the "Quality assurance" window.

Start of video instructions

As soon as one of the quality assurance or calibration procedures has been started, i.e. one of the corresponding buttons has been actuated, the related video instructions automatically start.

The steps to be performed are displayed in the respective video instructions as animated graphics with text, see Page 166.

The software internally checks whether all of the required inspection instruments are connected and in position. If this is not the case, the procedure just started cannot be continued.

Equipment delivery

Quality assurance

Treatment prescription

Treatment review

Calibration


Performing an isotropy check

- Actuate the <START> button.
The CC 600 executes an automatic procedure which takes several minutes:
 - X-ray radiation is generated and the CC 600 starts with the automatic procedure for checking the isotropy of the X-ray probe.
 - The five diagram bars of the displayed bar graph "Signal photo diodes [V]" (3) correspond to the five photo diodes in the PDA. The diagram bars are labeled +X, -X, +Y, -Y and +Z and display the current voltage on the respective photo diodes.
 - On expiry of the automatic procedure, the four bars, +X, -X, +Y, -Y, indicate the respective voltage for the optimal isotropy setting and the tolerance bars.
Since the tolerance bars are shown in reference to +Z, the height of the +Z bar remains unchanged with a deviation of 0V. If the X-ray radiation field is isotropic, the percentage deviation of the first four bars is near zero.
The tolerance bars on each of the diagram bars, +X, -X, +Y, -Y, indicate whether or not the isotropy is within the tolerance limits.

Legend for "Isotropy check" window

- 1 Button for <Isotropy check> function call
- 2 Graphic representation of radial beam position
- 3 "Signal photo diodes" bar graph, signals of five sensors
- 4 Text display of the current steps of the automatic process
- 5 Graphic representation of IRM counting rate (counts/time)
- 6 <START> button
- 7 <STOP> button

Result analysis

Isotropy out of tolerance

If the tolerance limits are exceeded, the color code "red" is displayed in the <Isotropy check> button. I.e., you must center the X-ray probe again. Change to the "XRS tip adjustment" function, Page 168.

Isotropy within tolerance

The XRS isotropy is within tolerance if the height of all five diagram bars is within the white tolerance bars and "green" is displayed in the color code field.

The isotropy reports can be displayed with the help of the "QA Reports" menu. These also are available as PDF files, see Page 181.

- Check the graphic representation (2) of the radial beam position (at the bottom left in the window):

The graphic representation displays the decentration of the electron beam in the X and Y directions. The red circle represents the absolute limit, i.e. if the decentration exceeds this value, the system switches off.


In practice, you should utilize maximally 50% of the tolerance, i.e. at roughly 50% of the maximum permissible decentration you should recenter the X-ray probe. To do this, call the "XRS tip adjustment" function, see Page 168.

The proper procedure for removing the X-ray probe from the inspection component (PAICH / PDA) is described in the section "Proper handling of components during quality assurance", see Page 165.

Latest test report (with PDA and PAICH)

<QA Reports> / <Open Latest>

Once you have carried out the quality assurance tests, you can display the results of the isotropy check (PDA) or output check (PAICH) that was last performed, i.e. the QA report that was most recently saved in your software.

Display test report, <Open Latest>

In the <QA Reports> menu, click <Open Latest> and select <PDA> or <PAICH>. The corresponding report is displayed as a table in a separate window.

Report browser

INTRABEAM 600

QA Isotropy Check

Institution:	Carl Zeiss Meditec AG
Date Time:	7.15.2016 11:26:05
Responsible Operator:	Service
XRS Serial Number:	507530
Beam Voltage [kV]:	50.0
Beam Current [uA]:	40.0
X-Offsets [mm]:	Static: 0.342 Dynamic: 0.420
Y-Offsets [mm]:	Static: -0.024 Dynamic: 0.260
PDA Model:	PDA
PDA Serial Number:	507571
Z Diode Output [V]:	4.617
+X Deviation [%]:	0.634
+Y Deviation [%]:	0.279
-X Deviation [%]:	1.244
-Y Deviation [%]:	0.181
Control Console Serial Number:	354
Result:	Pass

Done

Equipment inventory

Quality assurance

Treatment prescriptions

Treatment session

Calibration

Output check of dose rate

You must always perform this check!

With the "Output check" procedure, you determine the constancy of the dose rate of the XRS 4 using the IC ionization chamber and the dosimeter. In addition, the IRM rate is determined and compared with the calibration data.

The software calculates a correction factor between the measured dose rate and the data in the calibration file. This correction factor is automatically considered during the following treatment planning.

Successful implementation of the "Isotropy check" is a prerequisite for the "Output check". When the "Isotropy check" has been successfully completed, the <Isotropy check> button turns "green" and the <Output check> button turns "yellow".

Procedure

Probe adjuster	Recommended
Dynamic offsets	Recommended
Isotropy check	Mandatory
Output check	Mandatory

Function call

- Actuate the <Output check> button in the "Quality assurance" window.

Start of video instructions

As soon as one of the quality assurance or calibration procedures has been started, i.e. one of the corresponding buttons has been actuated, the related video instructions automatically start.

The steps to be performed are displayed in the respective video instructions as animated graphics with text, see Page 166.

The software internally checks whether all of the required inspection instruments are connected and in position. If this is not the case, the procedure just started cannot be continued.

Equipment inventory

Quality assurance

Treatment prescription

Treatment review

Calibration

Test setup for the output check

The proper procedure for inserting the X-ray probe in the PAICH using the V-guide and then setting up the XRS 4 in the XRS mount with the PAICH mounted vertically is described in the section "Proper handling of components during quality assurance", see Page 164.

Insert the X-ray probe in the PAICH using the V-guide.



WARNING

Danger due to defective ionization chamber

If you touch the sensor window, you may damage the ionization chamber up to the point of malfunction. This could possibly lead to overexposure.

- Do not touch the sensor window, touch the ionization chamber only on its side or on its cable.

The test setup is used to measure the dose rate currently being generated by the XRS 4 in the X-ray probe in the Z axis with the ionization chamber and the dosimeter. For the test setup, you must

- insert the XRS 4 in the PAICH and place it upright while positioning the XRS 4 exactly vertical in the XRS mount to prevent it from falling over,
- insert the ionization chamber in the mount on the PAICH and
- connect the cable of the ionization chamber to the "Ion Chamber" connector of the workplace. Make cable connection for XRS 4 and PAICH.

The dosimeter is automatically set to the measuring mode "Current", "Low pA" by the software. If you would like to change the dosimeter setting at this point, you must quit the software in order to access the dosimeter again for manual changes.

"Output check" procedure

The "Dosimeter status" color code display (5) is "green" when the dosimeter is communicating with the computer, i.e. when the dosimeter is active.

You must reset the dosimeter before beginning the "output check":

- Actuate the <Zero dosimeter> button (4).
- Actuate the <START> button (9) to start the procedure.
The procedure begins with a warm-up phase; then the values are read out and averaged.
- Check whether the ionization chamber is properly connected to the dosimeter, i.e. whether the ionization chamber is properly connected to the "Ion Chamber" connecting socket on the workplace.


I.a. the measured values of the dosimeter are displayed in the "Results" display field (7):

- "Dosimeter measurement [pA]", pure measurement
- "Dose rate output [Gy/min]", the value calculated internally taking the temperature and air pressure into account. The numerical value displayed is a time-averaged measured value of the dosimeter
- IBM rate and comparison with IRM calibration value.
- Actuate the <Accept> button (12) to confirm the measured value.

When the measured value has been confirmed, you can begin with the treatment planning (see Page 190) or perform a further quality assurance step.

Legend for "XRS output check" window

- 1 <Output check button
- 2 "QA completed" color code display
- 3 Air pressure display and <Unit> menu button for selecting unit of measure
- 4 <Zero dosimeter> button
- 5 "Dosimeter status" color code display

Equipment inventory

Quality assurance

Treatment prescription

Treatment review

Calibration

- Check the XY diagram on the bottom left of the screen mask which indicates the electron beam runout in the X and Y directions. The red circle represents the absolute limit value, i.e. the system is shut down when this value is exceeded.
In practical use, no more than 50% of this tolerance should be utilized. Therefore, you should recenter the XRS probe when approx. 50% of the admissible maximum runout are reached. Activate the "Probe Adjuster" function for this purpose, see Page 168.

The proper procedure for removing the X-ray probe from the inspection component (PAICH / PDA) is described in the section "Proper handling of components during quality assurance", see Page 165.

Equipment inventory

Quality assurance

Treatment prescription

Treatment review

Calibration


Output Check History, trailing window for "XRS output check"

The diagram shows the time variance of the XRS dose rate [Gy/min] and the percentage deviation referenced to the factory calibration over time [days since date of calibration file].

The "XRS Output % Deviation" display field (3) shows the percent deviation of the current measurement from the initial value (day 0).

If you actuate the <Yes> button (1), the measured values are saved.

If you actuate the <No> button (2), the measured values are rejected, i.e. you must perform the "XRS output check" once again.


Latest test report (with PDA and PAICH)

<QA Reports> / <Open Latest>

Once you have performed the "quality assurance" checks, you can display the respective results of the latest isotopy check (PDA) and output check (PAICH), i.e. the latest QA report most recently stored in your software.

Display test report, <Open Latest>

In the <QA Reports> menu, click <Open Latest> and select <PDA> or <PAICH>. The corresponding report is displayed as a table in a separate window.

INTRABEAM 600	
QA Output Check	
Institution:	Carl Zeiss Meditec AG
Date Time:	7.6.2016 11:17:22
Responsible Operator:	Ph. J. J. J.
XRS Serial Number:	207230
Beam Voltage [kV]:	50.0
Beam Current [nA]:	40.0
Dosimeter Model:	U4400-E
Dosimeter Serial Number:	000250
Dosimeter Calibration Factor:	1.00000E+0
Ion Chamber Model:	TN2342A
Ion Chamber Serial Number:	002691
Ion Chamber Calibration Factor [Gy.C]:	1.05000E+0
PAICH Model:	PAICH
PAICH Serial Number:	307372
Pressure [kPa]:	96.1
Standard Pressure [kPa]:	101.3
Temperature [°C]:	27.2
Standard Temperature [°C]:	22.0
P&T Correction Factor:	1.073
Output Value [Gy.min]:	5.892
Dose Rate Correction [a]:	1.0
IPM Count Rate [Hz]:	10000.495
IPM Change [a]:	-0.644
Control Console Serial Number:	354
Result:	Pass

Done

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page. There is no handwriting or other markings on the paper.

Equipment inventory

Quality assurance

Treatment prescription

Treatment approval

Treatment advice

Calibration

Treatment prescription

The "Treatment prescription" window is used to enter the treatment parameters for a specific treatment on a particular person.

The requirements for the "treatment prescription" are:

- The patient is registered in the patient database.
- The required components are registered in the equipment inventory.
- The quality assurance procedures were successfully completed.

Function call

- Actuate the <Treatment prescription> tab.

Selecting a patient

WARNING

Danger of confusion during treatment planning

If the user confuses the patient and/or the treatment plan, this could lead to irradiation being performed with the wrong treatment data. This could result in damage or injury due to overexposure.

- Make sure that you have selected the right patient and treatment plan.

Since you have to assign the treatment prescription to a patient, selection of the patient is the first operating step.

1 <Select patient> button (1).

When you actuate the <Select patient> button, you go to "Patient administration" and the "All patients" window is displayed, see Page 109:

- Select a patient, i.e. actuate the corresponding line in the table.
- Actuate the <OK> button, you then return to the "Treatment planning" window. The personal data of the selected patient are now displayed under "Patient information".

2 "Continue with previous patient" button


If the treatment of the selected patient is not performed for a period of more than 90 minutes, this button is "active".


- In order to continue treating the same patient within 90 minutes after the first treatment without a renewed inspection: Actuate the button (checkbox).

3 <Mark site> button

After you have selected a patient, define the treatment site in this treatment prescription:

- Actuate the <Mark site> button (2).
The "Treatment site" window is then displayed.
- For more information, see Page 192.

Marking the treatment site

Defining the treatment site

After the <Mark site> button is actuated in the "Treatment prescription" window, the "Treatment planning" window is displayed.

1 <Treatment site> menu button

Various different treatment sites are offered for your selection in a menu. The selected treatment site is displayed as a sketch. A cross hair which you can use to mark the treatment site is assigned to the sketch.

2 "Treatment site description" input field

You can define the treatment site by entering a description here.

3 <Apply> button

When this button is actuated, the current settings from the "Treatment planning" window are applied in the treatment prescription.

4 <Lock> button (lock icon)


When the lock icon is actuated, the open lock changes to a closed lock. This means that the data are entered and can be changed only after entering a password. The data can be edited, they are not archive data.


The data of a diagnosis are blocked; however, other treatment sites of the same patient are not affected.

5 X-axis

You can mark the position of the tumor with the cross hair. The cross hair can be moved manually.



6 Y-axis

You can mark the position of the tumor with the cross hair. The cross hair can be moved manually.

7 Keypad text: <Continue>

When this button is actuated, the current settings from the "Treatment site" window are applied in the treatment prescription.

Entering treatment parameters

The following requirements must be met for a treatment prescription:

- Selection of an applicator
- Selection of an XRS (X-ray source)
- Prescription of a dose

Fractional treatment

1 <Plan> button with status display

<Plan> color code for status display (plan for fractional treatment)

Green A current, complete plan for fractional treatment and the required equipment are already available.

Red The current plan is either incomplete or cannot be executed at present, because e.g. the equipment is inactive.

None Color code display:

No current plan for fractional treatment is available.

All previous fractional treatments have been completed

- Actuate the <Plan> button. The "Treatment planning" window is displayed. For more information, see Page 202.

Selection of an applicator

Only applicators for 50 kV are currently available. If you have prepared a treatment plan, this has already been taken into account there.


Only applicators which have been set to "Active" in the equipment inventory can be selected.

If the user selected an applicator for which no calibration data are available, the irradiation would be performed using incorrect treatment data.

2 <Type> menu button

The selection of a specific applicator type determines the type of treatment.

When you actuate the menu button all of the applicator types available are offered for selection in the component inventory. Select the required applicator type from the menu.

3 <Size> menu button


When you actuate the <Size> menu button, all sizes [cm] of the selected applicator type are offered for selection. Select the required applicator type from the menu.

4 <S/N> menu button

If several applicators of the same size are stored in the equipment inventory, you will need to select the required applicator according to its serial number.

If you tap the menu button <S/N>, the serial numbers of all the available applicators of the selected size [in cm] will be displayed for your selection.

- Compare the serial number of the applicator with the one in the software to ensure that you have selected the right applicator.

Equipment inventory

Quality assurance

Equipment approval

Treatment approval

Equipment status

Calibration

Selecting an XRS (X-ray source)

5 <S/N> menu button

When you actuate the <S/N> menu button, the serial numbers of all X-ray sources (XRS 4s) available in the equipment inventory are offered to you for selection. Only X-ray sources (XRS 4s) which have been set to "Active" in the equipment inventory and have an active quality assurance are offered for selection.

- Select an XRS 4 from the menu.
- Compare the serial number of the XRS with the one in the software to ensure that you have selected the right XRS.

6 <Voltage> menu button

The acceleration voltages at which the selected XRS 4 is calibrated are set out in a menu, from which you can select the desired voltage.

- Select the voltage you want to use to perform the treatment.

7 "Calibration mode" display

This displays the calibration mode that was used to calibrate the selected XRS 4 for the selected applicator.

For spherical applicators, you can choose between the calibration modes "TARGIT" and Calibration v4.0".

8 "Name" display



With the equipment inventory, you can assign a name to each XRS 4 in your INTRABEAM 600. Once you have selected an XRS 4 via its serial number, its name is displayed here.


9 <Current> menu button

The currents for which the selected XRS 4 has been calibrated are displayed in a table, from which you can select the desired current. Only currents which generate a minimum counting rate during the determination of the linearity of the IRM internal X-ray monitor (see Page 262) are entered in the calibration file and can be used for treatment planning.

- Select the current you want to use to perform the treatment.

Prescription (of a dose)

10 "Depth from surface" or "Depth to isocenter" input field (depending on the applicator selected)


- Enter the desired penetration depth (distance from the applicator surface or from the isocenter to the target area).

The numerical value entered can be edited using the "+" and "-" buttons located to the right of the input field.

11 "Dose" input field

- Enter the desired dose in the target area.
The dose in the target area results from the distance from the isocenter or from the distance from the applicator surface to the target area.

The numerical value entered can be edited using the "+" and "-" buttons located to the right of the input field.

12 "Treatment time" display field

When planning is completed, the software calculates the treatment time from the treatment parameters entered and displays them here.

13 "Dose rate at depth" display field

When planning is completed, the software calculates the depth dose rate from the treatment parameters entered and displays it here.


14 <Lock> button (lock icon)

The lock icon indicates whether the selected treatment plan is blocked (lock closed) or is not blocked (lock open).

The closed lock means that this treatment is protected for this patient and the treatment data cannot be edited in this window. The treatment data can be edited only after entering a password.

15 "Prescription complete" status display

Once you have entered all of the parameters for this treatment, this status display turns "green", i.e. the treatment planning is completed.

Equipment inventory

Quality assurance

Treatment prescription

Treatment approval

Treatment review

Calibration

Information on depth dose curve diagram

Once all parameters for a treatment have been entered in the "Treatment prescription" window, the depth dose curve (3) is graphically displayed in the lower part of the window.

The depth dose curve shows the dose [Gy] in relation to the distance from the applicator surface [mm] or to the isocenter [mm]. I.e. based on the depth dose curve you can determine what dose [Gy] occurs in the patient's body at certain depths [mm] (distances from the applicator surface).


Whether the distance of the applicator surface or of the isocenter is indicated depends on the applicator. With the needle applicator, it's the distance to the isocenter, otherwise it's the distance to the applicator surface.

It thus can be checked to what extent neighboring tissue is irradiated.

Two cross hairs are displayed in the diagram:

- The first, white cross hair (1) is fixed at the point of the treatment parameters entered during the prescription, i.e. the corresponding numerical values (coordinates) are entered in fields (6) and (7).
- The second, blue cross hair (2) can be moved manually. The current coordinates of the blue cross hair are displayed in fields (4) and (5) in the upper right corner of the diagram.

The "Depth" field (4) is a display and input field. You can enter an exact value for the depth [mm] here. The blue cross hair (2) is then positioned to this value in the diagram.

Planning fractional treatment

The treatment prescription for a patient also can be planned as a fractional treatment, i.e. the overall prescription can be distributed between several treatment appointments.

When you actuate the <Plan> button in the "Treatment prescription" window, the "Fractional treatment" window is displayed:

1 <Fractions> button

All of the treatment fractions saved for the current patient are displayed as a "Scheduled treatments" list.

The date of treatment is entered in the "Treatment date" column of the table for completed, administered fractions.

2 Selecting a fraction

- Actuate the respective line in the "Scheduled treatments" list. The selected fraction is color highlighted.

3 <Depth-Dose-Curve> button

The depth dose curve for the currently selected fraction is displayed in a separate window; for more information, see Page 200.

4 <Overview history> button

All changes to the fractional treatments are displayed in a separate window.

5 <Lock> button (lock icon)


The lock icon indicates whether the selected treatment plan is blocked (lock closed) or is not blocked (lock open).

The closed lock means that this treatment is protected for this patient. The treatment data can be edited only after entering a password.

6 <Add new fraction> button

You can create a new fraction (add a new line to the list). For a new treatment appointment, you must define an applicator, an XRS 4 and a dose prescription; for more information, see Page 194.


7 <Delete fraction> button

The currently selected fraction (color highlighted line in the "Scheduled treatments" list) can be deleted provided that it has not already been administered.

8 <Exit> button

Actuating this button returns you to the "Treatment prescription" window.

Blank lined paper with horizontal ruling lines.

Treatment approval

The "Treatment approval" is used to have the treatment parameters entered by the physicist verified and confirmed by the radiation therapist (treating physician). The radiation therapist (treating physician) confirms the correctness of the entered treatment parameters by entering his password.

"Treatment approval" window

The treatment parameters entered by the physicist via the "Treatment prescription" window are only displayed here and must be reviewed and confirmed by the radiation therapist (treating physician). Data can no longer be edited in the "Treatment approval". If data still have to be edited, you must return to the "Treatment prescription" (actuate the <Treatment prescription> tab).

The "Treatment approval" window can be called only if the "Treatment prescription" has been completed and the "Prescription complete" color code display is "green".

Function call

- Actuate the <Treatment approval> tab;
the "Treatment approval" window opens.

Password entry

Entering your password is the same as signing your name. Depending on the local legal situation, it may be necessary to have both the radiation therapist (treating physician) and the physicist sign the treatment record personally. The INTRABEAM 600 outputs the treatment record at the end of the treatment.

Password regulation

- If the physicist logs in to the INTRABEAM 600, this already is regarded as a signature.
- In case the physician logged in first, the following applies:
The physicist must sign and the physician must sign once again.

Only if the "Treatment approval" is completed correctly is it possible to actually perform the treatment, i.e. call the "Treatment monitoring" window.

"Treatment approval" window

1 Display of all treatment parameters

All treatment parameters (data) entered by the physicist in the current treatment prescription are only displayed here in the "Treatment approval" window and must be reviewed and confirmed by the radiation therapist (treating physician).

2 <Prescribing physician> menu button

This menu button is used to select the radiotherapist (treating physician) who is to check the current treatment prescription.

When the <Prescribing physician> menu button is actuated, the names of the physicians registered in your INTERBEAM 600 are offered for selection in a menu. Select the name of the physician by actuating the corresponding menu line.

3 "Password" input field

By entering his password, the radiation therapist (treating physician) confirms the treatment approval of the currently displayed treatment prescription.

4 <Verify> button

By actuating the <Verify> button, the radiation therapist (treating physician) confirms that the entered treatment parameters are correct, thus concluding the treatment approval of the current treatment prescription.

5 "Verification complete" color code display

Once the radiation therapist (treating physician) has actuated the <Verify> button, the color code "green" is displayed by the "Verification complete" indicator light, i.e. the inspection of the treatment parameters for the current treatment prescription is concluded. This is a prerequisite for subsequently being able to call the "Treatment monitoring" window.

245 | Maier, Jörg

1975-04-05 | Male | 2.000 kg | brain

Fracture | Depth-Dose-Curve | Parameters | Mark treatment site

	Site	Physician	Site Locked
1	2016-01-27	Physician	Site Locked
2	2016-01-27	Physician	not
3	2016-01-27	Physician	Site Unlocked
4	2016-01-27	Physician	Site Locked

Exit

Operation

Checklist for the INTRABEAM 600	210
Conducting treatment	212
"Treatment monitoring" window	216
Starting the treatment.....	218
Interrupting the treatment (Pause mode)	222
Canceling a treatment (Stop)	224
Treatment terminated.....	226
Troubleshooting.....	227
 Review treatment, displaying treatment parameters	231
Displaying the treatment record.....	242
Generating a treatment record	244

Checklist for the INTRABEAM 600

- Always check the following points before each use of the device (without the patient!):

The "OR equipment checklist" can be found in the chapter "Supplementary documents", Page 316.

Operational readiness

- ✓ The calibrations of the components required for the current treatment are all still valid.
- ✓ The workplace is fully equipped with the components and accessories required for the current treatment.
- ✓ The INTRABEAM 600, XRS 4 and stand have been verified and confirmed.
- ✓ All of the components required for the current treatment have been properly cleaned and disinfected.

Accessories

- ✓ The additional equipment, for example INTRABEAM Applicators, INTRABEAM Drape etc., have been checked based on their respective instructions for use, are sterile and are available as complete sets.
- ✓ The expiry date has not lapsed.
- ✓ The calibration data of the applicators have been loaded into the application software.

External interlock

- ✓ Door contact for external interlock connected and ready for operation
- ✓ "X-ray" warning lamp connected and ready for operation

Device readiness

- ✓ The workplace is properly positioned, i.e. it is located far from the operating table (patient) and outside of the sterile operating area.
- ✓ The workplace has been secured against rolling away via the locking devices on the front wheels.

- ✓ The workplace and the stand are connected to the power supply.
- ✓ The workplace and the stand are switched on.
- ✓ A cable connection has been made from the XRS 4 to the workplace.
- ✓ The stand is properly positioned, balanced and secured against rolling away by locking the brakes of the stand base.
- ✓ XRS 4 with applicator properly seated in the mount of the stand.
- ✓ The INTRABEAM Drape is attached.
- ✓ The stand has been balanced correctly.
- ✓ The XRS 4 with mounted applicator are correctly positioned in the patient according to the treatment prescription.

Equipment inventory

Quality assurance

Treatment prescription

Treatment planning

Treatment review

Calibration

Conducting treatment

WARNING

Danger due to unshielded radiation

If the generation of X-rays is started before the X-ray probe has been brought into the target area, unshielded radiation will be emitted, thus harming the patient and the user.

- Do not switch the X-ray source on until the X-ray probe with applicator are located in the target area.

WARNING

Radiation hazard for critical organs

If critical organs are not adequately shielded during treatment, they may be damaged by X-ray radiation.

- Use radiation shields to shield critical organs.

WARNING

Danger due to radiation emerging from the patient.

If the patient does not offer sufficient shielding against emerging radiation or sufficient arrangements for radiation protection have not been made on the patient, X-ray radiation may emerge from the patient into his surroundings and harm the user.

- Provide adequate shielding in the target area using the radiation shields.
- If necessary, protect personnel with portable radiation shields.
- Keep as far away from the target area as possible.

WARNING

Danger due to radiation during treatment

During irradiation, X-ray radiation can emerge from the patient into his surroundings, thus harming the user and other persons in the treatment room.

- It is essential for you to equip your treatment room with an external interlock (interlocking device) and an X-ray display (warning lamp).
- Before beginning irradiation, ensure that all persons who are not wearing suitable radiation protection suits have left the treatment room.

Equipment inventory

Quality assurance

Treatment prescription

Treatment monitoring

Treatment review

Calibration

**CAUTION**

The probe tip of the XRS X-ray source contains incompatible material. If the user brings the X-ray probe directly into contact with the patient's tissue, incompatible material contacts the wound tissue, thus possibly causing allergic reactions.

- Never insert an X-ray probe directly into the patient's tissue.
- Always place an applicator on the X-ray probe.

Tumor bed irradiation

- Switch on the stand and the workplace.
- Mount the XRS 4 with probe protector on the stand.
- Connect the cable from the XRS 4 mounted on the stand to the workplace. (see page 103)
- Mount the applicator on the stand determined by the treatment prescription.
- Always cover the stand sterily with an INTRABEAM Drape. The proper procedure for covering the stand sterily with an INTRABEAM Drape is described in the Instructions for Use of the "NC32 INTRABEAM Floor Stand".

**CAUTION****Danger due to cross-contamination**

Failure to use a drape can cause cross-contamination in the sterile operating area.

- Always use the stand in the operating room only with an INTRABEAM Drape.

NOTE**Heat build-up during treatment**

One inevitable side effect of generating X-ray radiation is the simultaneous generation of thermal energy in the X-ray probe of the XRS 4 and the dissipation of thermal energy into the patient's surrounding tissue.

- Please note that: The surface temperature of applicators may exceed 41°C during the course of treatment.

Equipment inventory

Quality assurance

Treatment preparation

Treatment monitoring

Treatment review

Calibration

Calling the "Treatment monitoring" window

You perform irradiation using the "Treatment monitoring" window.

Before calling the "Treatment monitoring" window, you must ensure that your INTRABEAM 600 is completely set up; in particular, all connector and cable connections must have been made properly.

The "Treatment monitoring" window can be called only if the "Treatment approval" has been completed, i.e. the "Verification complete" color code display in the "Treatment approval" window is "green", see Page 205.

- Ensure that the stand with the XRS 4 and mounted applicator is positioned correctly according to the currently pending treatment.

Function call

- Actuate the <Treatment monitoring> tab;
The "Treatment monitoring" window is displayed.

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page. There is no handwriting or other markings on the paper.

Equipment monitoring

Quality assurance

Treatment prescription

Treatment monitoring

Treatment record

Calibration

"Treatment monitoring" window

1 Patient information

The treatment parameters from the treatment prescription are listed here. They comprise information on:

- the individual patient
- the diagnosis
- the applicator
- the XRS 4

2 Running treatment

You will find the most important information on the running or current irradiation treatment here:

- Dose planned
- Dose delivered, i.e. dose already emitted
- IRC counts (IRM counting rate), percentage specification as diagram
- Beginning, start time of irradiation
- Remaining time
- Elapsed treatment time, percentage indication as diagram

3 "Treatment status" display field

The "Treatment status" display field reads "Preparation" until the treatment parameters have been loaded to the CC 600. Then the status display changes to "Ready".

4 <Treat> button

The button is initially called <Treat>. If the <Treat> button is actuated, the <Treatment intro> window opens and, when you exit this window, the text of the button changes to <Load data>. When this button is actuated, the treatment parameters are loaded to the CC 600. When the loading process has been completed, the status display changes to "Ready" and the button changes to <Start>. The irradiation starts.

5 <STOP> button, for canceling the treatment

6 Interlock status

The interlock status indicates whether all interlock contacts of the INTRABEAM 600 are closed and what state they are in.

7 Display and comment field

Depending on your selection, either the event log or an input field for comments is displayed here.


8 <Comments> button

When you actuate the <Comments> button, the comment field is shown in the display and comment field (7). You can enter and edit your comments on the current treatment here.

Your comments are automatically saved. After the treatment, the comments can be displayed in the "Treatment review" window. Then your comments no longer can be edited there.

9 <Log File> button

When you actuate the <Log File> button, the event log is displayed in the display and comment field (7). The event log is a display field.

Starting the treatment

The INTRABEAM 600 checks all interlock contacts and all cable connections. If all conditions are not fulfilled, the button (4) is inactive, i.e. the status display for this button is set to "red".

- If all conditions for the treatment are fulfilled, i.e. the status displays are "green" and the button (4) is set to <Ready>, the treatment can be started:
- Actuate the button (4) in the <Ready> setting.

Checking components

Prior to irradiation, all of the prescribed components are displayed in a separate window along with their serial number and the interlock status.


WARNING

Danger of confusion when selecting components

If a component whose data does not match the data entered in the treatment description is used for the treatment, this can lead to radiation damage or injuries.

- Compare the reference number, the serial number and the UDI (unique device identification) to make sure that the data and the component are correctly assigned.

Copyright notice

Quality assurance

Copyright notice

Copyright notice

Copyright notice

If the serial numbers on the current components match the ones from the treatment prescription, you may actuate the <Continue> button.

As soon as the CC 600 is ready, the designation of the button in the "Treatment monitoring" window (4) changes from <Treat> to <Load data>.

Loading treatment data to the CC 600:

It must be ensured that the treatment data are saved in case of a communication failure between the computer and the CC 600 (e.g. due to a computer crash).

- Actuate the <Load data> button.
After a delay of approx. 2 seconds, the treatment data are loaded on the CC 600 and the word <Start> appears on the keypad.

Now you must not make any more changes to your INTRABEAM 600; e.g. you must not undo any plug connections.

- Actuate the <Start> button (4) to begin the irradiation.

At the start of irradiation, the status displays change to "XRS on" and the designation of the button (4) changes from <Start> to <Pause>.

NOTE

Heat build-up during treatment

One inevitable side effect of generating X-ray radiation is the simultaneous generation of thermal energy in the X-ray probe of the XRS 4 and the dissipation of thermal energy into the patient's surrounding tissue.

- Please note that: The surface temperature of applicators may exceed 41°C during the course of treatment.

Graphic representation of course of treatment


10 Information on diagram (legend)

11 Diagram (curve progression)

The standardized IRM rate over the running time of the treatment is shown as a curve with marked upper and lower limits.

12 Graphic representation of radial beam position

Graphic representation and indication of numerical value of current beam deflection and indication of current X and Y position.

Treatment history

Quality assurance

Treatment parameters

System status

Calibration

- **Monitor the following data during the treatment:**

The current values of the IRM internal X-ray monitor should be within a tolerance of $\pm 10\%$ of the specified values. If the current values are outside this tolerance, the application software reports an error and cancels the treatment.

- Observe the beam deflection in the "Radial beam position" field (12): If the current values are out of tolerance, the irradiation is interrupted and an error message appears.

The current beam position usually drifts slowly to the circle perimeter. The reason for this is often tension of the X-ray probe caused by a change in the position of the patient or the tissue. A short pause to reduce the mechanical stress of the X-ray probe can prevent a cancellation of the treatment.

4 <Pause> button

During irradiation, the <Pause> button is used to interrupt the treatment. During the pause, the word <Continue> appears on the keypad and is used to resume the treatment.

The button is initially called <Treat>. If the <Treat> button is actuated, the <Treatment intro> window opens and, when you exit this window, the text of the button changes to <Load data>. When this button is actuated, the treatment parameters are loaded to the CC 600. When the loading process has been completed, the status display changes to "Ready" and the button changes to <Start>. The irradiation starts.

Interrupting the treatment (Pause mode)

You may require a treatment pause to perform a minor intervention on the patient, e.g. to correct his or her position or if a large deflection is detected in the graphic representation of the radial beam position (12). During the treatment pause you also can, for example, easily correct the position of the XR 4 in order to reduce the force applied to the X-ray probe.

A brief pause to relieve mechanical tension on the x-ray probe can avoid the need to interrupt the treatment.

- Actuate the <Pause> button (4) to interrupt (pause) the current irradiation. A long beep sounds.

"Pause" is displayed in the Status display field. During the pause, the button (4) is designated <Continue>.

Since the INTRABEAM 600 does not emit any radiation during the pause, the curve of the IRM rate (11) must drop zero.

You can actuate the <Continue> button (4) to continue the treatment. When you continue the treatment, the designation of the button (4) again changes to <Pause>.


NOTE

Treatment cannot be interrupted

If a treatment cannot be interrupted via software because, for example, the <Pause> button does not react when it is actuated, this may result in an impairment of the patient and/or damage to the XRS 4.

- Actuate the <PAUSE> button on the control panel of the CC 600, see Page 70.

This button also can be used to interrupt (pause) the treatment. However, on the CC 600, repeated actuation of the <PAUSE> button can not be used to continue the treatment.

Equipment inventory

Quality assurance

Treatment preparation

Treatment execution

Treatment review

Calibration


Canceling a treatment (Stop)

Treatment can be canceled: By actuating the <STOP> button (5) or the Emergency Off button on the control panel of the CC 600 (see "Control panel of the CC 600" Page 70).

Where possible it is preferable to cancel the treatment by using the <STOP> button (5), because this means it is still possible to save all of the treatment data in the proper manner. If the Emergency Off button is used, for technical reasons only the most important data can be saved.

The INTRABEAM 600 displays a corresponding window indicating the dose emitted so far and its deviation from the prescribed nominal dose in percentage. (The numerical values shown here are only examples.)

WARNING

Treatment cannot be cancelled

If a treatment cannot be canceled via software control because, for example, the <STOP> button (5) does not react when it is actuated, this may lead to overexposure.

Actuate the Emergency Off button on the control panel of the CC 600, see Page 70.

In order to switch the CC 600 back on, you first must unlock the Emergency Off button (by turning it and pulling it out).

After canceling a treatment, you must perform another quality assurance check of the XRS 4 (isotropy check and output check) unless you continue treating the previous patient within 90 minutes.

Equipment inventory

Quality assurance

Treatment prescription

Treatment monitoring

Treatment review

Calibration

Continue with previous patient

In order to be able to continue treating the same patient within 90 minutes, you have to go to the "Treatment prescription" of the current treatment, see Page 190.

- Actuate the <Treatment prescription> tab.
- Actuate the "Continue with previous patient" button.
In the "Treatment prescription" window the data of the current treatment are displayed if the treatment of the current patient has not been canceled for more than 90 minutes.


Treatment terminated

The treatment is terminated according to plan when the delivered irradiation dose is equal to the prescribed irradiation dose (see Treatment planning). If the delivered irradiation dose deviates from the prescribed irradiation dose, for example in case of a cancellation, the software provides support with the window, see Page 224. A fraction with the remaining dose can be generated by the software automatically. Also applies to "Troubleshooting".

Continuing treatment of the same patient

If you would like to treat the same patient again, you can plan an additional treatment within 90 minutes without having to verify the XRS 4 again.

- To do this, place a check mark next to "Continue with previous patient" in the "Treatment prescription" window, see Page 225.

Troubleshooting

NOTE

Following a failure, the user may perhaps be unable to power up the INTRABEAM 600 again without any information from the Instructions for Use and then have to cancel the treatment.

- Ensure that the users receive special training regarding troubleshooting and the Instructions for Use are always within reach.

Failure of the computer

If the data exchange between the CC 600 and the computer is interrupted while the treatment is in progress, e.g. due to a crash of the computer software or for similar reasons, the CC 600 then continues the treatment alone.

If the user cancels the treatment due to failure of the computer, then the length of the treatment will be unnecessarily extended.

- Let the INTRABEAM 600 continue working; do not in any case restart the INTRABEAM software.

If the CC 600 continues the treatment alone, all of the information required for the treatment will be output on the control panel display of the CC 600:

Control panel of the CC 600 during the treatment:

	Nominal value (setting)		Nominal irradiation dose [Gy]		
Acceleration voltage [kV]	50 kV	S	5.00 Gy	9.50 m	Nominal irradiation period [m], in minutes
Current intensity [µA]	40 µA	A	2.50 Gy	4.75 m	Actual irradiation period [m], in minutes
	Actual value		Actual irradiation period		

On reaching the set target dose, the CC 600 automatically terminates the treatment.

- When the CC 600 has terminated the treatment, restart the computer. When the computer restarts, the CC 600 automatically exchanges treatment data with the computer, so you can check and output all of the data from the "Treatment review" functional group on the computer.

Failure of the CC 600

If the CC 600 is disconnected from the power supply during treatment or switched off by actuation of the Emergency Off button, the treatment data generated up to that point are automatically stored. Then proceed as follows:

Shut down the software properly with <File>/ Exit

Switch the CC 600 on again.

Wait until the self-test has been completed and restart the software.

The CC 600 and the computer then automatically exchange treatment data so that you can check and output all important data in the "Treatment review" functional group of the software on the computer.

In this case, less than the planned dose was emitted. The differential dose can be delivered with new specifications in a further treatment. If the differential dose is to be delivered with the INTRABEAM 600, re-inspection of the XRS 4 is required or, if desired, irradiation can be continued without reverification using the "Continue with previous patient" function.

In cases where the CC 600 must be replaced, the computer creates a "preliminary" treatment record. This record is replaced by a final record later on following renewed communication with the CC600.

Equipment delivery

Quality assurance

Treatment prescription

Treatment monitoring

Treatment review

Calibration

Automatic pause mode

If an interlock error occurs during the treatment, for example:

- The applicator no longer is properly seated in its mount on the stand,
- An external interlock (door contact) is interrupted,
- An XRS cable has been removed on the XRS 4 etc. ...,

the INTRABEAM 600 will then automatically go into the Pause mode. In the "Treatment monitoring" window "Pause" is displayed as the current operating status.

- Close all interlocks completely until the "Interlock status" display for "External" and "XRS" again turns "green".
- After this, press <Continue>, on the keypad to continue with the treatment.

Troubleshooting for the ABC error

(ABC error = automatic beam centering error)

An error message may occur if the probe tip of the XRS 4 is deflected too strongly. The reason for this is often tension of the X-ray probe caused by slight changes in the position of the patient or the tissue.

If the current values of the beam position are out of tolerance, irradiation is automatically interrupted and an error message appears.

If mechanical tensions act on the applicator during treatment, you can observe this in the "Treatment monitoring" window on the graphic representation of the radial beam position and determine that the beam position slowly drifts towards the circle perimeter.

- In this case, stop the treatment with the <Pause> button.

The mechanical tension acting on the applicator can be eliminated by briefly releasing the brakes of the stand.

NOTE

Danger when brakes of stand are released

Releasing the brakes of the stand may result in an uncontrolled movement (i.e. lowering) of the stand arm. This can injure the patient and damage the X-ray probe.

A sterile force is required to hold the stand arm in place while the brakes are briefly released.

If an error occurs and the software does not offer automatic generation of a fraction with the residual dose:

- Note the previously delivered and the prescribed dose.
- Determine the differential dose.
- Use the hand switch to release the brakes on the stand in order to reduce the stress of the probe tip of the XRS 4.
- Actuate the "Treatment prescription" tab, plan this differential dose with "Continue with previous patient" as a new treatment and continue irradiation.

If the error occurs again:

The XRS 4 must be removed from the treatment field by a sterile person.

- Actuate the "Quality assurance" tab.
- Realign the X-ray probe, i.e. carry out the XRS tip adjustment.
- Perform all other quality assurance routines, i.e. the dynamic offset, isotropy check and output check.
- Plan a new treatment with the calculated remaining dose.

If the error persists:

- Inform ZEISS Service.

Miscellaneous treatment discontinuation

- Note the previously delivered and the prescribed dose.
- Determine the differential dose.
- Actuate the "Treatment prescription" tab, plan this differential dose with "Continue with previous patient" as a new treatment and continue irradiation.

Symptoms summary

Quality assurance

Treatment parameters

Evaluation

Review treatment, displaying treatment parameters

Using the "Treatment review" window, the course of a treatment can be traced and the treatment record can be created.

The parameters of the currently terminated treatment are automatically displayed in the "Treatment review" window. Records of older treatments also can be called up.

Treatment record of the currently completed treatment

When a treatment is called up, the "Treatment review" automatically displays a summary of the data from the currently completed treatment.

These data are automatically saved in a database where the respective file name is automatically generated from the "patient information". For more information, see Page 120.

Function call

- Open the "Treatment review" screen mask by clicking the <Treatment review> tab.

Treatment record of a previous treatment

Saved records of previous treatments also can be displayed in the "Treatment review".

- Actuate the <Load> button.
The <Record selector> window is displayed.

Multiple treatments sorted by treatment ID are listed for selection in the <Record selector> window.

- Select a treatment by actuating the corresponding line.
- Actuate the <Select> button.
The detailed treatment data are displayed in the <Treatment review> window.

If a treatment is loaded into the <Treatment review> window, a treatment record can be created and saved using various different invocations, i.e. <Single treatment>, <Short form> or <Multiple fractions>.

Equipment inventory

Quality assurance

Treatment form updates

Treatment review

Calibration

The procedure for having the treatment records exported is described in the section "Editing folders, exporting files" on Page 122.

Information on the "Treatment review" window

The treatment parameters of the currently loaded treatment are displayed in the upper part of the "Treatment review" window (Items 1 to 6).

1 Patient information

The information of the patient being treated is displayed.

2 Applicator

The applicator currently being used is displayed.

3 XRS

The XRS 4 being used is displayed.

4 Responsible team

Radiotherapy physicist and physician who are responsible for this treatment

5 Treatment duration

"Duration
(delivered)"

Display of the actual treatment duration of the irradiation performed. The actual treatment time is without any possible pause times. If the dose rate changes during the treatment, the actual treatment time may deviate from the planned treatment time by maximally 10%.

"Duration
(planned)"

Display of the originally planned treatment duration for comparison.

6 Prescription

The planned and the actually delivered dose as well as the penetration depth of this treatment prescription are displayed.

7 Display field

Various different types of information from the treatment record can be displayed using the tabs of this display field.

8 <Multiple fractions> button

Selection of a fractional treatment and creation of a corresponding treatment record, for more information see Page 244.

9 <Short form> button

Generate short form of treatment record from the treatment currently displayed in the "Treatment review" window. For more information, see Page 248.


10 <Single treatment> button.

Generate treatment record of a <Single treatment> from the treatment currently displayed in the "Treatment review" window. For more information, see Page 250.

11 <Load> button

A saved treatment record of a completed single treatment can be loaded for the "Treatment review". When you actuate the <Load> button, the program branches to the "Load record" window. Here you can select the desired single treatment by actuating the <Select record> button; for more information, see Page 242.

Display field with tab bar

The display field (7) of the "Treatment review" window displays various information about the treatment that has just been loaded. Select the desired information type by actuating the corresponding tab in the tab bar above the display field.

12 Tab of <Comments> display field

The comments entered manually during the treatment are shown in the display field.

13 Tab of <Beam position> display field

The chronological sequence of the beam deflection is shown in the display field.

14 Tab of <IRM> display field

The current conditions (expected / actual) of the IRM internal X-ray monitor counting rate are shown in the display field.

15 Tab of <Depth dose> display field

In the display field, the depth dose curve (depth [mm] / dose [Gy]) of the planned and the actually provided dose is displayed.

16 Tab of <Log file> display field

In the display field, the record of the treatment is displayed as a list, i.e. as an event log (currently shown in the display field).

Equipment inventory	Quality assurance	Treatment prescription	Treatment review	Calibration
---------------------	-------------------	------------------------	------------------	-------------

16	15	14	13	12	7
File Name: 1017-1 Depth dose: 1017-1 BM: 1017-1 Beam position: 1017-1 Comments: 1017-1	Equipment inventory Equipment: 1017-1 Equipment ID: 1017-1 Equipment Name: 1017-1 Equipment Type: 1017-1 Equipment Status: 1017-1 Equipment Location: 1017-1 Equipment Date: 1017-1 Equipment Time: 1017-1 Equipment User: 1017-1 Equipment Password: 1017-1 Equipment Comment: 1017-1	Quality assurance Quality ID: 1017-1 Quality Name: 1017-1 Quality Type: 1017-1 Quality Status: 1017-1 Quality Location: 1017-1 Quality Date: 1017-1 Quality Time: 1017-1 Quality User: 1017-1 Quality Password: 1017-1 Quality Comment: 1017-1	Treatment prescription Treatment ID: 1017-1 Treatment Name: 1017-1 Treatment Type: 1017-1 Treatment Status: 1017-1 Treatment Location: 1017-1 Treatment Date: 1017-1 Treatment Time: 1017-1 Treatment User: 1017-1 Treatment Password: 1017-1 Treatment Comment: 1017-1	Treatment review Treatment ID: 1017-1 Treatment Name: 1017-1 Treatment Type: 1017-1 Treatment Status: 1017-1 Treatment Location: 1017-1 Treatment Date: 1017-1 Treatment Time: 1017-1 Treatment User: 1017-1 Treatment Password: 1017-1 Treatment Comment: 1017-1	Calibration Calibration ID: 1017-1 Calibration Name: 1017-1 Calibration Type: 1017-1 Calibration Status: 1017-1 Calibration Location: 1017-1 Calibration Date: 1017-1 Calibration Time: 1017-1 Calibration User: 1017-1 Calibration Password: 1017-1 Calibration Comment: 1017-1

27-Jan-16
14-16-20 Treatment screen selected
Responsible Physician/Physician
Responsible Physician/Physician
Patient ID: 245
Treatment No: 1017-1

14:19:03 Patient check passed
14:19:03 READIR state reported by the Control Console
14:19:03 All Interlock conditions are satisfied
14:19:37 Loading Treatment Parameters into the Control Console, please wait ...
14:19:40 Treatment Parameters loaded
14:19:41 READIR state reported by the Control Console
The treatment is ready to start.

View treatment report
Generate treatment report


1-1017
Load
Single treatment
Short form
Multiple fractions

Displaying comments

The manually entered comments on the currently selected treatment can be displayed in the display field (7).

- Actuate the tab of the <Comments> display field (12).


Displaying the beam position

The beam position (offset in X and Y direction [mm] over the running time [min:sec]) of the currently selected treatment can be displayed in the display field (7).


- Actuate the tab of the <Beam position> display field (13).

Displaying the IRM counting rate

The IRM counting rate of the currently selected treatment can be displayed in the display field (7). This is the curve progression of the standardized counting rate of the internal X-ray monitor [%] over the running time [min:sec].

- Actuate the tab of the <IRM> display field (14).


Displaying the depth dose curve

The depth dose curve of the currently selected treatment can be displayed in the display field (7). The specified depth dose curve from the treatment prescription and the depth dose curve actually provided are compared.

- Actuate the tab of the <Depth dose> display field (15).

Displaying the event log

The event log of the currently selected treatment can be displayed in the display field (7).

- Actuate the tab of the <Log file> display field (16).

The screenshot shows the INTRABEAM 600 interface with the following components:

- Navigation Bar:** Equipment inventory, Quality assurance, Treatment prescription, **Treatment review** (selected), Calibration.
- Display Field (7):** Shows the event log for the selected treatment. The log includes the following entries:
 - 14:10:03 Patient check passed
 - 14:10:03 READY state reported by the Control Console
 - 14:10:03 All interlock conditions are satisfied
 - 14:10:37 Loading Treatment Parameters into the Control Console, please wait
 - 14:10:40 Treatment Parameters loaded
 - 14:10:42 READY state reported by the Control Console
- Log file tab (16):** The selected tab for displaying the event log.
- Other tabs:** Equipment inventory, Quality assurance, Treatment prescription, Calibration.
- Buttons:** View treatment report, Generate treatment report, Load, Single treatment, Short term, Multiple fractions.

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page. There is no handwriting or other markings on the paper.

Displaying the treatment record

You can display a saved treatment record of a completed single treatment in the "Treatment review" window. For this, you must load this treatment record in the "Treatment review" window so that it can be displayed there. Use the <Load> button to call the trailing window, "Record selector". In this window you will find a list of all the treatment records stored in your INTRABEAM 600. You can select the desired treatment record from this list.

Function call

- Actuate the <Load> button in the "Treatment review" window. The "Record selector" window is then displayed.

"Record selector" window

- Select a single treatment (line in the list) by actuating the relevant line. The line is color highlighted.
- Actuate the <Select> button.
The "Record selector" window closes and the record of the selected single treatment is displayed in the "Treatment review" window.

<Cancel> button

Without having selected a single treatment, close the "Record selector" window by actuating the <Cancel> button.

[illegible]

Record number	Record ID	Record Name	Record Date	Record Time
1-1052	123456	Elisa Mustekau	2014-03-14	12:29:56
1-1054	123456	Elisa Mustekau	2014-03-13	09:40:48
1-1055	123567	Alan Mustekau	2014-03-15	09:50:33

Generating a treatment record

Generating a treatment record of a fractional treatment

Based on the data from fractional treatments saved in your INTRABEAM 600, you can generate the respective treatment records using the "Record selector" function.

Use the <Multiple fractions> button to call the "Record generator" window. This window contains two lists. You can select a patient in the upper list. The fractional treatments previously performed for the selected patient are displayed in the lower list. Select the fractional treatment from which you would like to generate a treatment record.

Function call

- Actuate the <Multiple fractions> button (8) in the "Treatment review" window. The "Record generator" window is then displayed.

"Record generator" window

- Select a patient in the upper list (line in the list) and a "treatment site" if several sites are listed.
Actuate the relevant line for this purpose. The line is color highlighted.
- Select a treatment fraction in the lower list (line in the list) by actuating the relevant line. The line is color highlighted.

<  > button (sum function)

Using this button, you can simultaneously select multiple treatment fractions and generate a joint treatment record from them.

If a treatment fraction (line in the list) has been selected and you then actuate the <  > button and subsequently select additional treatment fractions, all of them are jointly selected, i.e. the selected lines are all color highlighted.

The selection can be canceled (deactivation of the "sum function") by actuating the <  > button again.

<Generate PDF> button

A treatment record is generated in PDF format from the selected treatment fractions and displayed.

You can save it in the "Report browser" when then opens.

To save the report, please only use the <Save PDF> button.

<Cancel> button

Without having generated a single treatment record, close the "Record generator" window.

Record generator

Record generator

Record ID	Patient	Diagnosis	Date of birth	Sex	Age	Center
123567	Max	Melanoma	02-06-1976	Male	070543	2222

Show all

Fr	Seq	Beam	Beam ID	Beam Type	Beam Energy	Beam Shape	Beam Size	Beam Date	Beam Time	Beam Physician	Beam Physician	
1	1-1053	1.00	1.00	0.0	0.0158	507530	A61005	Spherical	10-12-98	09-14-2016	Physician	Physician
2	1-1054	1.00	0.93	44.7	0.0104	507126	A61005	Spherical	08-08-98	08-15-2016	Physician	Physician

Generate PDF Cancel

Equipment inventory

Legal & assurance

Treatment prescription

Treatment group

Location

Report browser


INTRABEAM 600

Fractional Treatment Report

Institution: Hospital Oberkochen
 Report prepared by: Service
 Report generation date: 15-Jul-16 13:04:39
 Patient Given Name: Max
 Patient Family Name: Mustermann
 Patient Date of Birth: 1960-02-02
 Patient Sex: Male
 Patient ID: 1
 Diagnosis: glioblastoma
 Insurance ID:
 Other ID:
 Overall prescription dose [Gy]: 40.000
 Overall delivered dose [Gy]: 40.000
 Overall deviation from prescription [%]: 0.000
 Prescription Depth [mm]: 0.00

#	Trt No.	Plan [Gy]	Actual [Gy]	Delta [%]	Trt Time	NRS S/N	APP S/N	APP Type
1	354-66	40.00	40.00	0.0	0:41:06	507530	F61065	0.0 Spherical
#	Trt No.	Started	Date	Physicist	Physician			
1	354-66	11:42:10	07-06-2016	Physicist	Physician			

Save PDF

Done

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page. There is no handwriting or other markings on the paper.

Treatment identity

Quality assurance

Treatment preparation

Treatment review

Calibration

Generating a short form treatment record

Based on the data of the currently concluded treatment, you can generate a short form treatment record using the "Short form" function.

Function call

- Actuate the <Short form> button (9) in the "Treatment review" window. In a separate window, the short form treatment record of the currently loaded treatment is displayed.

Control functions of the "Record generator" window:

1 <Done> button

Without having saved the currently displayed treatment record, terminate the "Short form" function; the "Treatment review" window then reappears.

2 <Export reports> button

The currently displayed treatment record can be exported to DICOM.

3 <Save PDF> button

The currently displayed treatment record can be saved as a PDF file in your INTRABEAM 600. You should only use the <Save PDF> button to save the report; otherwise, it cannot be made secure.

Equipment inventory

Quality assurance

Treatment prescription

Treatment review

Calibration

Report browser

INTRABEAM 600

Treatment Record

Institution	aaa
Patient Given Name	Max
Patient Family Name	Mustermann
Patient Date of Birth	10-Okt-2000
Patient Sex	Male
Patient ID	123456
Diagnosis	Brain
Treatment No.	1-1045
Insurance ID	Abc321
Other ID	123abC
Treatment Date	01-Mrz-18
Treatment Start/Stop Time	11:58:58 / 12:02:05
Prescribing Physician	Physician
Responsible Physicist	Physicist
Control Console SN	1
XRS Serial Number	507528
XRS Name	2nd
Calibration Mode	Target
Voltage [kV]	50.0
Current [µA]	40.0
Applicator Type	Spherical
Applicator Size [cm]	5.0
Applicator Serial Number	H01065
Applicator Sterilization: to Date	4
Applicator Sterilization: Max	160
Estimated Treatment Time	0:02:05
Actual Treatment Time	0:02:06
Prescription Dose [Gy]	0.50
Prescription Depth [mm]	5.00
Dose Delivered [Gy]	0.50
Deviation from Prescription [%]	0.0
FDA Report Name	OutputCheck_XRS507528@50.0kV_20160302114554.pdf
PAICH Report Name	IsotropyCheck_XRS507528@50.0kV_20160302113638.pdf

Seite 1 von 18

Save PDF

Export report

Done

☐ Anonymous report

4

3

2

1

Generating a treatment record of a single treatment

Based on the data of the currently concluded treatment, you can generate a treatment record using the "Single treatment" function.

Function call

- Actuate the <Single treatment> button (10) in the "Treatment review" window. In the "Record generator" window, the treatment record of the currently concluded treatment is displayed.

Control functions of the "Record generator" window:

1 "Anonymized report" checkbox

The currently displayed treatment record can be separated from the personal information by actuating the checkbox.

A check mark is displayed in the checkbox and the personal information is no longer displayed in the treatment record. Actuating the checkbox again resets the anonymization.

Anonymized treatment records cannot be exported to DICOM, i.e. the <Export reports> button remains inactive as long as the treatment record is anonymized.

2 <Done> button

Without having saved the currently displayed treatment record, terminate the "Single treatment" function; the "Treatment review" window then reappears.

3 <Export reports> button

The currently displayed treatment record can be exported to DICOM.

4 <Save PDF> button

The currently displayed treatment record can be saved as a PDF file in your INTRABEAM 600. You should only use the <Save PDF> button to save the report; otherwise, it cannot be made secure.

Equipment inventory

Quality assurance

Treatment prescription

Treatment record

Calibration

Report browser

INTRABEAM 600

Treatment Record

Institution	333
Patient Given Name	Max
Patient Family Name	Mustermann
Patient Date of Birth	10-Okt-2000
Patient Sex	Male
Patient ID	123456
Diagnosis	Brain
Treatment No.	1-1045
Insurance ID	Abc321
Other ID	123abC
Treatment Date	02-Mrz-16
Treatment Start/Stop Time	11:59:58 / 12:02:05
Prescribing Physician	Physician
Responsible Physician	Physician
Control Console SN	1
XRS Serial Number	507528
XRS Name	2nd
Calibration Mode	Target
Voltage [kV]	50.0
Current [nA]	40.0
Applicator Type	Spherical
Applicator Size [cm]	3.0
Applicator Serial Number	H01065
Applicator Sterilization: to Date	4
Applicator Sterilization: Max	100
Estimated Treatment Time	0:02:05
Actual Treatment Time	0:02:06
Prescription Dose [Gy]	0.50
Prescription Depth [mm]	5.00
Dose Delivered [Gy]	0.50
Deviation from Prescription [%]	0.0
PDA Report Name	OutputCheck_XR5507528@50.0kV_20160302114554.pdf
PAICH Report Name	IsotropyCheck_XR5507528@50.0kV_20160302113638.pdf

Seite 1 von 18

Save PDF

Export report

Done

☐ Anonymous report

4

3

2

1

Blank lined paper with horizontal ruling lines.

Care and maintenance

Calibration instructions	254
Calibration	257
Determining the beam centering.....	258
Determining the linearity of the IRM internal X-ray monitor	262
Calibrating the air pressure sensor.....	264
Calibrating the temperature sensor	267
Safety check.....	270
Workplace safety check.....	270
Troubleshooting table.....	271
Changing the line fuse of the workplace	277
Handling, transport and shipment.....	278
Handling components of the INTRABEAM 600	278
Care of the INTRABEAM 600	279
Handling components	281
Cleaning components	282
Disinfecting components	284
Maintenance intervals.....	285
Disposal	288

Equipment inventory

Quality assurance

Treatment preparation

Treatment review

Calibration

Calibration instructions

The following calibrations are required for safe and reliable operation of your INTRABEAM 600:

Calibrating the XRS 4

During a calibration, the characteristic of the depth dose rate of the X-ray field is quantified for the respective XRS 4, i.e. the output and the radiation characteristic are calibrated.

The calibration interval begins with the initial startup. Since several months may possibly pass between the calibration and the initial startup, biannual calibration must be set on initial delivery. Nevertheless, ZEISS will automatically call your attention to the expiry of the calibration interval within one year following initial startup.

- Reset the calibration interval to "Annual" when the first calibration interval expires.

WARNING

Danger due to noncompliance with calibration deadlines

If the calibration deadlines of the XRS 4 are not met, this can lead to incorrect irradiation.

- Have your XRS 4 calibrated annually (as recommended by ZEISS). If national regulations require shorter calibration intervals, these must be observed.
- Have your XRS 4 calibrated by ZEISS at the manufacturer's plant. As part of this calibration, a new calibration file is created for the XRS 4 and delivered together with the XRS 4.

Logistical matters

Quality assurance

Regulatory compliance

Document control

Calibrating the dosimeter and the ionization chamber

Use of the ionization chamber (IC) on the dosimeter (UNIDOS E from the PTW company) requires regular calibration of the ionization chamber and the dosimeter.

The PTW company recommends calibration every two years.

If national regulations require shorter calibration intervals, these must be observed.

WARNING

Danger due to measuring errors of the ionization chamber and the dosimeter

If the calibration intervals of the ionization chamber and the dosimeter are not observed, the quality assurance measurements may be incorrect. This could possibly lead to under- or overexposure.

- Have your ionization chambers and your dosimeters calibrated at least every two years as recommended by the PTW company.

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There is no text or other markings on the paper.

Calibration

Initial calibrations are performed by the manufacturer, ZEISS, at the manufacturer's plant. Generally, the INTRABEAM 600, its individual components and its accessories are always calibrated on delivery.

You can perform additional calibrations yourself with the help of the software.

ZEISS offers calibrations at the manufacturer's plant. Observe the calibration regulations, Page 254.

"Calibration" function call

- Click the <Calibration> tab.
The "Calibration" window is displayed.

"Calibration" window

You can call up the following functions in the "Calibration" window:

- Procedure
- 1 —•
 - 2 —•
 - 3 —•
 - 4 —•

1 Beam centering

Determine the sensitivity of the beam deflection of an X-ray probe. see Page 254.

2 IRM linearity

The linear correlations between the individual beam currents and the corresponding IRM counting rate in relation to the counting rate for the maximum beam current are determined for the IRM internal X-ray monitor. see Page 262.

3 Barometer

Compare the air pressure measured by the CC 600 with a calibrated external barometer and correct if necessary. see Page 264.

4 Thermometer

Compare the temperature measured by the PAICH with a calibrated external thermometer and correct if necessary. see Page 267.

Equipment inventory

Quality assurance

Treatment prescription

Treatment review

Calibration

Determining the beam centering

You can use this function to perform a static offset calibration, i.e. to determine the sensitivity of the deflection coils of an XRS 4. For this purpose, you can record a characteristic curve which shows the voltage change that is caused in the probe tip by a certain change in the beam direction of the electron beam.

For the arrangement of the components and the procedure for adjusting the "Dynamic offset" beam deflection, see Page 173.

You can import the newly determined sensitivity of the deflection coils shown below in the "Deflection factor curve" window to the calibration file of the XRS X-ray source concerned.


We recommend performing this calibration at least once annually and on receipt of a new calibration file.

"Calibration" function call

- Click the <Calibration> tab.
The "Calibration" window is displayed.
- Select <Beam centering> under "Procedure".

Procedure

Start of video instructions

As soon as one of the quality assurance or calibration procedures has been started, i.e. one of the corresponding buttons has been actuated, the related video instructions automatically start.

The steps to be performed are displayed in the respective video instructions as animated graphics with text, see Page 166.

The software internally checks whether all of the required inspection instruments are connected and in position. If this is not the case, the procedure just started cannot be continued.

Equipment inventory

Quality assurance

Treatment prescription

Treatment review

Calibration

Operation of "Calibrate beam deflection"


<START> button

Actuate the <START> button to start an automatic measuring procedure.

The measuring procedure takes several minutes. At the end of the automatic measuring procedure, the measured values are displayed in the "Deflection Factor History" window and automatically imported to the calibration file of the XRS 4 concerned.

<STOP> button

You can cancel the automatic measuring procedure with the <STOP> button.


Equipment inventory

Quality assurance

Equipment prescription

Treatment review

Beam centering history

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page. There is no handwriting or other markings on the paper.

Equipment inventory

Quality assurance

Treatment prescription

Treatment setup

Determining the linearity of the IRM internal X-ray monitor

"Calibration" function call

- Click the <Calibration> tab.
The "Calibration" window is displayed.
- Select <IRM linearity> under "Procedure".

Procedure

IRM linearity

The linear correlations between the individual beam currents and the corresponding IRM counting rate in relation to the counting rate for the maximum beam current are determined for the IRM internal X-ray monitor. I.e. you can determine the linear dependency of the IRM counting rate, starting from the maximum current (40 μA). With an exactly linear behavior, the counting rate at 10 μA would be exactly 25%.

This step is not required for an acceleration voltage of 40 kV, since treatment at 40 kV is possible only with **40 μA** .

We recommend performing this calibration at least once annually and on receipt of a new calibration file.

Only currents which generate a minimum counting rate can be used for treatment planning.

Start of video instructions

As soon as one of the quality assurance or calibration procedures has been started, i.e. one of the corresponding buttons has been actuated, the related video instructions automatically start.

The steps to be performed are displayed in the respective video instructions as animated graphics with text, see Page 166.

The software internally checks whether all of the required inspection instruments are connected and in position. If this is not the case, the procedure just started cannot be continued.

Using the function "Determining the linearity of the IRM"


<START> button

Actuate the <START> button to start an automatic measuring procedure.

The measuring procedure takes several minutes. At the end of the automatic measuring procedure, the measured values are displayed in the "IRM linearity" window and automatically imported to the database for the respective XRS 4.

<STOP> button

You can cancel the automatic measuring procedure with the <STOP> button.

Equipment inventory

Quality assurance

Treatment prescriptions

Treatment review

Calibration

Calibrating the air pressure sensor

The INTRABEAM 600 requires the exact air pressure for the correction calculation of the dose during the "output check of the dose rate", see display field in the window "Output check of XRS dose rate", Page 182.

An air pressure sensor is integrated in the CC 600. In order to calibrate this air pressure sensor, you need a calibrated external barometer. During calibration, the CC 600 and the calibrated external barometer must be located in the same room. Define the calibration interval according to the locally applicable statutory provisions.

WARNING

Danger due to lack of or improper calibration of the air pressure sensor

Failing to observe the calibration interval or setting the pressure units [mm Hg, kPa or bar] inconsistently in the software or on the calibrated barometer can result in an incorrectly calculated dose rate. This could possibly lead to overexposure.

- Calibrate your air pressure sensor at least once a year.
- Ensure that the air pressure values in the software and on the calibrated barometer have the same units of measurement.

"Calibration" function call

- Click the <Calibration> tab.
The "Calibration" window is displayed.
- Select <Barometer> under "Procedure".

Operation of the "Barometer calibration" window


1 "Internal reading" display

This displays the air pressure value measured by the air pressure sensor integrated in the CC 600.

2 "Method" menu button

As a user of the INTRABEAM 600, you are entitled to perform a "1 Point Calibration".

Procedure

Equipment inventory

Quality assurance

Treatment prescription

Treatment review

3 "Unit" menu button

Make sure that the pressure unit you set on your INTRABEAM 600 is the same as that in which the air pressure is displayed on the external barometer.

4 "Enter external reading" input field


Read the current air pressure reading on the calibrated external barometer and enter this numerical value here. The software calculates a new compensation factor.

5 <Cancel> button

With <Cancel> you exit the window without saving a new compensation factor.

6 <Save> button

The newly calculated compensation factor is saved with the <Save> button.

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page. There is no handwriting or other markings on the paper.

Equipment inventory


Quality assurance

Treatment prescription

Treatment review

Calibration

Calibrating the temperature sensor

WARNING

The INTRABEAM 600 requires the exact temperature for the correction calculation of the dose during the "output check of the dose rate", see display field in the "Output check" window, Page 182. If an incorrect temperature is determined, the dose correction calculation results in an incorrect dose value. Follow the instructions for calibrating the temperature sensor.

In the PAICH; a temperature sensor is integrated in the mount for the IC ionization chamber. In order to calibrate this air pressure sensor, you need a calibrated external thermometer.

- Place the temperature sensor of the calibrated external thermometer in the PAICH, in the mount for the IC ionization chamber.
- Wait until the temperature reading displayed by the calibrated external thermometer is stabilized, i.e. until the temperature reading practically no longer changes.

Danger due to lack of or improper calibration of the temperature sensor


If the calibration intervals are not observed or differently set units of measurement [°C] or [°F] are used for temperature values in the software and on the calibrated thermometer, an incorrectly calculated dose rate will result. This could possibly lead to overexposure.

- Calibrate your temperature sensor at least once a year or define the calibration interval according to the locally applicable statutory specifications.
- Ensure that the temperature values in the software and on the calibrated thermometer have the same units of measurement.

"Calibration" function call

- Click the <Calibration> tab.
The "Calibration" window is displayed.
- Select <Thermometer> under "Procedure".

Procedure

Equipment inventory

Quality assurance

Feedback prescription

Treatment review

Calibration

Start of video instructions

As soon as one of the quality assurance or calibration procedures has been started, i.e. one of the corresponding buttons has been actuated, the related video instructions automatically start.

The steps to be performed are displayed in the respective video instructions as animated graphics with text, see Page 166.

The software internally checks whether all of the required inspection instruments are connected and in position. If this is not the case, the procedure just started cannot be continued.

Operation of the "Thermometer calibration" window

1 "Internal reading" display

The temperature reading of the temperature sensor integrated in the PAICH is displayed here.

2 "Method" menu button

As a user of the INTRABEAM 600, you are entitled to perform a "1 Point Calibration".

3 "Unit" menu button

Make sure that the temperature unit you set on your INTRABEAM 600 is the same as that in which the temperature is displayed on the external thermometer.

4 "Enter external reading" input field

Read the current temperature reading on the calibrated external thermometer and enter this numerical value here. The software calculates a new compensation factor.

5 <Cancel> button

With <Cancel> you exit the window without saving a new compensation factor.

6 <Save> button

The newly calculated compensation factor is saved with the <Save> button.

The diagram illustrates the process flow for the preparation of the test specimen, organized into four vertical columns (1, 2, 3, 4) representing different stages or equipment groups. The process steps are as follows:

- Step 1:** Preparation of the test specimen. This step involves the use of equipment from column 1 and column 2.
- Step 2:** Preparation of the test specimen. This step involves the use of equipment from column 2 and column 3.
- Step 3:** Preparation of the test specimen. This step involves the use of equipment from column 3 and column 4.
- Step 4:** Preparation of the test specimen. This step involves the use of equipment from column 4.

The diagram also includes a 'Flowchart' section with a 'Flowchart' title and a 'Flowchart' description. It includes a 'Flowchart' section with a 'Flowchart' title and a 'Flowchart' description. The diagram is a complex flowchart showing the steps involved in the preparation of the test specimen. It includes a 'Flowchart' section with a 'Flowchart' title and a 'Flowchart' description. The diagram is a complex flowchart showing the steps involved in the preparation of the test specimen. It includes a 'Flowchart' section with a 'Flowchart' title and a 'Flowchart' description.

Safety check

Workplace safety check

To prevent any impairment to the safety of the workplace as a result of aging, wear, etc., the organization operating the system must ensure, in accordance with the applicable national regulations, that the regular technical safety checks defined for this workplace are performed on schedule and to the stipulated extent.

The safety checks must only be performed by the manufacturer or qualified personnel.

At least the following safety checks should be performed for the workplace:

- Availability of the Instructions for Use
- Visual inspection of the workplace, the accessory components and the cable connections for damage and legibility of the labels
- Function test of the folding or drawer function of the respective compartments
- Leakage current test
- Test of protective ground conductor
- Function and wear test of double casters and their brakes
- Function test of all switches, buttons, sockets and indicator lights.

Troubleshooting table

Although the components of the INTRABEAM 600 are high-quality technical products, malfunctions and error messages may nevertheless occur. The following "Troubleshooting table" will support you if a malfunction or an error message occurs.

WARNING

The INTRABEAM 600 does not contain any components that the user himself can service. If an unauthorized person attempts to perform service on the components of the INTRABEAM 600, uncontrolled X-rays may be emitted and components may be damaged. This could possibly lead to overexposure.

- If a malfunction or error message occurs that you cannot correct according to the "Troubleshooting" table, attach a sign to the device stating that it is out of order and contact ZEISS Service.
- Do not open any components of the INTRABEAM 600.
- Have service work performed on your INTRABEAM 600 only by ZEISS Service or by a specialist designated and authorized by ZEISS.

Service Identification Program (SIP)

Your device is registered in the ZEISS "Service Identification Program" (SIP). Your device bears an SIP label with an SIP number, see Page 50. The SIP number is a unique identification number assigned to your device or system. ZEISS maintains a product file under this SIP number.

- Please send your service requirements to your local contact of the ZEISS sales branch in your country. The telephone number, fax number and e-mail address are all listed on the SIP label.

Malfunction or error message	Possible cause	Remedy	See
No function at all.	Power plug of stand not plugged in.	• Plug in the power plug.	
	Power failure.	• Contact in-house electrician.	
	Power switch of workplace not switched on.	• Set power switch to "I". Green indicator light in power switch must light up.	
	Fuses defective.	• Disconnect workplace from power grid and replace fuses in the power supply module.	
The computer doesn't start when the workplace is switched on.	Malfunction during autostart of computer.	• Switch workplace off, wait for 30 seconds and switch it on again.	
CC600 display is empty.	CC600 not functioning.	• Call ZEISS Service.	
Communication between CC 600 and computer interrupted.	One or more cables are loose.	• Call ZEISS Service.	
	Computer did not boot up properly.	• Switch workplace off, wait for 30 seconds and switch it on again.	Page 111
Software is no longer operable		• Restart the computer.	Page 112
		• In case of repeated recurrence, contact ZEISS Service.	
Keyboard and/or mouse not working.	USB cable of keyboard and/or mouse defective.	• Call ZEISS Service.	

Malfunction or error message	Possible cause	Remedy	See
CC600 Error during self-test	The INTRABEAM stand and the workplace are not connected via the XRS cable – one end of the cable is not connected.	• Connect both ends of the XRS cable • with the INTRABEAM stand • with the workplace • Turn off the computer and wait until the screen goes black. • Switch off the workplace, wait for 30 seconds and then restart it.	
ID 10015 You must work OFFLINE, since the answering CC 600 is not listed for this purpose in the equipment inventory. The answering CC 600 outputs its serial number on the monitor." You can enter the CC 600 in the equipment inventory and then log in again under Main menu -> Tools -> Login.	The software has detected that the serial number of the currently connected CC 600 does not agree with the serial number of the CC 600 activated in the equipment inventory.	• Start with "Offline" in order to set the control unit currently being used to "Active" in the equipment inventory. • If the current CC 600 has not yet been entered and stored in the equipment inventory, then do so immediately and restart the software.	Page 116 Page 140
ID 10016 The dosimeter was not found by the software! Please make sure that the dosimeter is switched on and the serial cable is connected. Select 'Repeat' when you're finished or select 'Continue' to continue without a dosimeter. In this case, you can switch the dosimeter on later.	The RS232 communication cable between the PC and the dosimeter is not connected properly.	• Check dosimeter for proper connection.	Page 102 --

Malfunction or error message	Possible cause	Remedy	See
ID 10018 The answering dosimeter is not 'Active' in the equipment inventory. The answering dosimeter outputs its serial number of the monitor. Continue without a dosimeter or set the dosimeter to 'Active' in the equipment inventory. Then you can log in again via the Main menu -> Tools -> Login.	The software has detected that the serial number of the currently connected dosimeter does not match the serial number of the dosimeter active in the equipment inventory.	• Start "Offline" and set the dosimeter currently being used to "Active" in the equipment inventory. • If the current dosimeter has not yet been entered and stored in the equipment inventory, then do so and subsequently restart the software. • Exchange the dosimeter according to the displayed serial number and connect the RS232 cable. • Switch dosimeter on.	Page 116 / Page 148 Page 148 Page 148 Page 66
ID 001 Invalid user ID or invalid password! Please try once again or exit the program.	A false password and/or user ID was entered.	• Repeat login with correct user ID and password.	Page 114
ID 201 Please confirm the XRS 4 calibration data.	New XRS 4 calibration data have been imported.	• Click <OK> and check whether the serial number and data are correct. • Click <Accept> to accept the new data. • Click <Reject> to reject the new data.	-- -- --

Malfunction or error message	Possible cause	Remedy	See
ID 202 Please confirm the applicator calibration data.	New applicator calibration data have been imported.	- Click <OK> and check the serial number and the data. - Click <Accept> to accept the new data. - Click <Reject> to reject the new data.	-- -- --
ID 203 Previously available *.cal file missing.	Original *.cal file no longer located in the directory provided for that purpose.	Click <OK>. The file in question may be located in the archive directory, seeing as a newer .cal file has been loaded. Or the file was renamed or deleted.	--
ID 10019 Filament not ready. Please wait ...	The XRS 4 is not at operating temperature.	Wait until the XRS 4 reaches operating temperature. Start the procedure again if necessary.	--
ID 10021 The CC 600 is now performing irradiation. Please wait until the irradiation has been completed. Exit the program or go offline.	Computer is trying to communicate with a CC 600 in the treatment mode. An unintended constellation is involved here.	- Wait until the treatment is finished and do not use the computer. - The end of treatment is indicated on the control unit display and when the "X-Ray On" LED goes out. - Restart the software.	Page 227 -- --
IRM Rate Error 0x00000008	- Defect in the XRS cable or - Hardware error in CC 600 or - fault in XRS	Call ZEISS Service.	

Malfunction or error message	Possible cause	Remedy	See
"0x00000400 high voltage error" "0x00000800 beam current error" (Malfunction may typically occur shortly after the start)	- Defect in the XRS cable or - Hardware error in CC 600 or - fault in XRS	Call ZEISS Service.	--
"0x00100000 control error"	Irradiation was interrupted because the automatic electron beam centering no longer was able to compensate for the decentration of the tip of the XRS X-ray source.	- Check the position of the XRS X-ray source for mechanical stress and correct if necessary. - If you want to continue irradiation, plan a new irradiation for the still missing dose with the data in the "Treatment review" window. In this case, you must inspect the XRS 4 again. You also can use a different, previously prepared XRS 4.	--- Page 231 Page 158 ff Page 140
0x00001000 IRM PIC error (IRM PIC: Monitoring of the IRM rate on the mainboard of the CC 600). Expected IRM rate out of tolerance.	- Plug connectors on the XRS cable or - The XRS cable itself or - Hardware fault in the CC 600 or - Fault in the XRS 4.	- Check both plug connections of XRS cable. - Continue treatment by selecting the function: "Continue with previous patient".	Page 277 Page 191
0x00004000 VIT PIC error Control unit problem with evaluation of signals	Control unit hardware fault.	Call ZEISS Service.	

Changing the line fuse of the workplace


The line fuses are located in the power supply module (1) on the back of the workplace.

The order data of the replacement fuse can be found on page 292.

If necessary, you also can have the line fuse changed by our Service or by an authorized person.

The following procedure is mandatory when changing the line fuse:

- Switch the workplace on at the power switch.
- Unplug the power plug from the workplace.
- Use a screwdriver with a narrow tip to release the clips of the fuse holder (2).
- Pull out the fuse holder.
- Remove the defective line fuses from the fuse holder.
- Insert the new line fuses in the fuse holder.
- Slide the fuse holder into the fuse drawer until the clips re-engage (2).
- Reconnect the workplace to the power grid.
- Switch the workplace on at the power switch.

Handling, transport and shipment

Handling components of the INTRABEAM 600

- Observe the following when handling, transporting and shipping components of the INTRABEAM 600:
 - When handling contaminated components of the INTRABEAM 600, always wear protective clothing such as gloves, an apron or a surgical gown and eye protection.
 - Always assume that a component of the INTRABEAM 600 is contaminated if it is not known whether it is clean, disinfected or sterile.
 - For transport of the INTRABEAM 600 or its individual components, make sure that the parts being transported are properly protected against vibrations and strong temperature fluctuations.
 - Transport the components of the INTRABEAM 600 inside the clinic only in the transport and sterilization containers provided for this purpose or transport the complete INTRABEAM 600 (with control unit, computer and dosimeter) in the workplace.

NOTE**Risk of functional deterioration**

If components of the INTRABEAM 600 are transported or sent by post over long distances in unsuitable packaging, they may be damaged.

- Transport or send components of the INTRABEAM 600 outside of the clinic only cleaned and disinfected in their original packaging or in a shuttle container provided by ZEISS.

Care of the INTRABEAM 600

CAUTION

Danger due to improper processing

If components, applicators or parts not intended for processing are processed by means of processes not approved for this purpose, their asepsis or sterility cannot be guaranteed and damage may result. This could possibly lead to overexposure.

- Use only the methods approved for each respective component for processing.
- In order to protect patients and operating personnel against infections, use only components which have been carefully cleaned and disinfected during treatment.
- Use only sterile applicators.
- Always cover the stand with a sterile disposable cover (INTRABEAM Drape).
- Do not recycle any disposable products.
- Use only the approved procedures described here to process the components.

For further details and explanations on the approved cleaning and disinfection procedures, please refer to Page 282 ff.

The sterilization procedures for sterilizable applicators are provided in the separate, respective Instructions for Use of the applicator concerned.

The following INTRABEAM 600 components and accessories must be cleaned and disinfected in compliance with the respective care instructions:

- Workplace
- V-guide (workplace)
- Monitor (touchscreen monitor)
Observe the special instructions for the touchscreen!
- Keyboard (membrane keypad)
- Mouse
- XRS 4
- PDA
- PAICH

- XRS and PDA/PAICH cables
- Stand
- INTRABEAM Water Phantom
- Dosimeter*
- IC ionization chamber*

* For instructions on cleaning the UNIDOS E dosimeter and the IC ionization chamber, please refer to the Instructions for Use from the PTW company, Freiburg.

The following applicators can be steam sterilized:

- Spherical applicators and set with transport and storage containers
- Surface applicator set
- Flat applicator set

The sterile disposable products are:

- INTRABEAM Needle Applicator
- Radiation Shield, Flat 12"x17"
- INTRABEAM Drape

Handling components

When handling the components of the INTRABEAM 600, note the following:

NOTE**Danger in case of nickel allergy**

If a user with a nickel allergy has skin contact with the housing of an XRS 4 or an inspection instrument (PAICH, PDA), allergic reactions may result.

- Always wear suitable gloves and protective clothing if necessary.

**CAUTION**

The probe tip of the XRS 4 contains beryllium (Be). The outer surface of the probe tip is coated with a protective layer. If this protective layer is damaged, a person's skin and clothing may be contaminated.

- Avoid direct skin contact with the beryllium and send in the XRS 4 to the manufacturer immediately.
- First aid measure: Following skin contact with beryllium, wash with plenty of water. Remove contaminated clothing.

NOTE

If the components are not stored in the appropriate transport and storage containers, they may be damaged and contaminated, and thus rendered unusable.

- Always store the components in the appropriate transport and storage containers if they are not being used.

NOTE

Dirty or tarnished interlock surfaces can lead to situations in which the interlock does not perform its function and the INTRABEAM 600 cannot generate X-ray radiation.

- Clean the interlock surfaces only with the approved cleaning agents.

Cleaning components

WARNING

Danger due to illegible symbols or labels on the device!

Improper cleaning methods or wear may make symbols and/or labels on your device illegible or cause them to come undone. This can lead to reading errors and confusion. Use of the wrong components could possibly lead to overexposure.

- If you notice that any symbols or labels are missing on your device or have become illegible, please contact us or one of our authorized representatives.

We will supply you with replacements.

NOTE

Damage caused by penetrating moisture

Penetrating moisture damages electronic components.

- Make sure that no water or cleaning agent seeps into the electronic components:
 - Do not spray electronic components with cleaning or disinfectant solutions.
 - Wipe off electronic components only with a damp cloth.

The protective CrN coating of the X-ray probe can be damaged.

- Do not leave the XRS 4 x-ray probe sitting in water for longer than 2 hours. Failure to comply with this instruction may result in damage to the protective CrN coating of the X-ray probe.
- If possible, clean the mechanical surfaces of the components and accessories immediately after use.
- Do not allow contaminations should to dry on the objects, as this would make cleaning and disinfecting more difficult.
- Check the entire surfaces of the individual components for coarse contaminations and other foreign particles.
- Remove any contaminations or foreign substances by wiping them off with a moist disposable sponge or with gauze pads moistened with water.

- Do not use any aggressive or abrasive cleaning agents.
- Remove residue using a mixture of 50% ethyl alcohol and 50% distilled water plus a dash of household dish-washing liquid.
- Wipe off the surfaces moist repeatedly until all visible contaminations have been removed. You can use an abrasion-free cleaning cloth if necessary.
- Use a lint-free swab or a soft-bristle brush instead of a gauze pad when cleaning hard-to-access areas (in tubes or pipes).
- Let the components dry before continuing to use them.

Cleaning the touchscreen

- Please note:
 - Do not use any all-purpose cleaners, other cleaning agents or scouring agents.
 - No moisture may enter the touchscreen. Therefore, do not spray water or glass cleaner directly onto the touchscreen, but wipe it on with a cloth instead.
- Switch the INTRABEAM 600 off.
- Wipe off the monitor screen and the housing with a soft, clean and dry cloth (e.g. a microfiber cloth). If necessary, moisten a cloth slightly with water or a pure glass cleaner and wipe off the touchscreen with the moistened cloth.

Disinfecting components

It may be necessary to disinfect surfaces (especially those of control elements).

NOTE**Risk of damaging surfaces!**

- Use a disinfectant based on aldehyde and/or alcohol. The addition of quaternary compounds is acceptable. To prevent damaging surfaces, disinfecting components other than those listed below must not be used.

The maximum concentrations are:

For alcohol (tested with 2 propanol): 60%

- For aldehyde (tested with glutaraldehyde): 2%
 - For quaternary compounds (tested with DDAC): 0.2%
- Let the components dry before continuing to use them.

NOTE**Damage caused by penetrating moisture**

Penetrating moisture damages electronic components.

- Make sure that no water or cleaning agent seeps into the electronic components:
 - Do not spray electronic components with cleaning or disinfectant solutions.
 - Wipe off electronic components only with a damp cloth.

**The protective CrN coating of the X-ray probe can be damaged.**

- Do not leave the XRS 4 x-ray probe sitting in water for longer than 2 hours. Failure to comply with this instruction may result in damage to the protective CrN coating of the X-ray probe.

Maintenance intervals

Regularly checking and maintaining your device is necessary to ensure safe and correct operation of the system and a long service life. The following table lists the relevant maintenance activities and shows the maximum intervals at which they must be performed.

If the device has to be opened for maintenance work, this may be done only by authorized specialists. Please contact your local ZEISS Service in due time for the performance of these maintenance activities. You will find the contact authorized for your country on the following website:
www.meditec.zeiss.com

If you find out during self-performed tests that the device no longer meets its acceptance criteria, you must not operate the device any longer! Please shut down the device and contact your local ZEISS Service.

Regardless of the maintenance work specified, operators of medical electrical devices are obliged to perform periodic inspections in compliance with IEC 62353.

Component	Inspection features and acceptance criteria
No regular inspections or maintenance calls are required in order to maintain electro-magnetic compatibility (EMC).	
According to country-specific interval requirements	
Safety check	Check the device safety, radiation protection and electrical safety according to the prevailing national regulations.
Metrological control	According to national requirements
Recalibration of measuring devices	Recalibration of the dosimeter (UNIDOS E) and ionization chamber according to the manufacturer's specifications
Prior to each use (by the user)	
Stand (NC32 Floor Stand)	Electrical and mechanical function test: Release axis brakes completely; braking force holds XRS securely in balanced state without any change of position.

Component	Inspection features and acceptance criteria
XRS 4, PAICH, PDA, Ionization chamber	Visual inspection for external damage
Every 12 months (by authorized service personnel as per maintenance plan)	
Metrological inspection of XRS 4	Metrological inspection and maintenance by the manufacturer or authorized organizations
Documentation and identification labels	Visual inspection: Instructions for Use and identification labels available, undamaged, complete and legible
Controls, handgrips, buttons	Visual inspection and function test: controls operate smoothly and without play.
Stand (NC32 Floor Stand)	Function test: – Play-free movement of bearings, stops and suspensions – Axis brakes hold XRS 4 firmly and securely.
Casters	Visual inspection and function test: – Play-free seating, no locking during rotating and rolling – Secure hold of the locking devices and brakes – Cable guards present
Fans of CC 600 and workplace	Function test: Fan performance noticeable.
Touchscreen	Visual inspection and function test: – No noticeable damage to housing – Check image quality – Check functioning of joints
CC 600	Check of optical and acoustic signaling. Check of correct air pressure display
Every 2 years (by authorized service personnel as per maintenance plan)	
Electronics	Battery change – Computer (battery on mainboard) – CC 600 (if equipped with a battery)

Component	Inspection features and acceptance criteria
-----------	---

Every 4 years (by authorized service personnel as per maintenance plan)	
--	--

Stand (NC32 Floor Stand)	Visual inspection: – No signs of wear on the supporting structure – No signs of wear or damage on the interfaces and couplings
Cables, plugs and switches	Visual inspection: – Undamaged insulation and connectors on device/ cables – Check of cable routing including fastenings

Every 6 years (by authorized service personnel as per maintenance plan)	
--	--

Casters	We recommend replacing the casters

Every 7 years (by authorized service personnel as per maintenance plan)	
--	--

Electronics	Since the expected service life of electronic storage devices is approx. 7 years, it may be necessary to replace electronic components after this period.
-------------	---

Disposal

User information on the disposal of electrical and electronic devices

This symbol means that the product must not be disposed of as normal domestic waste.

The correct disposal of electrical or electronic devices helps to protect the environment and to prevent potential hazards to the environment and/or human health which may occur as a result of improper handling of the devices concerned.

For detailed information on the disposal of the product, please contact your local dealer or the device manufacturer or its legal successor. Please also note the manufacturer's topical information on the Internet. In the event of resale of the product or its components, the seller is required to inform the buyer that the product must be disposed of in accordance with the applicable national regulations currently in force.

For end customers in the European Union

Please contact your dealer or supplier if you wish to dispose of electrical or electronic devices.

Information on disposal in countries outside the European Union

This symbol is only applicable in the European Union. For the disposal of electrical and electronic devices, please observe the relevant national legislation and other regulations applicable in your country.

Disposal of outer packaging material

- Dispose of outer packaging material in accordance with local regulations.

System data

Order data	290
Accessories	291
Spare parts	292
Technical data	293
Technical data of the INTRABEAM Workplace	293
Technical data of isolating transformer	294
Technical data of the XRS 4	295
Technical data for UNIDOS E	296
Technical data for Control Console 600 (CC 600)	297
Technical data of computer	297
Ethernet connection (optional),	298
Data carriers must be virus-free.	299
External interlock and X-radiation indicator	300
External X-ray display	302
EMC (electromagnetic compatibility)	303
EMC - Electromagnetic compatibility as per IEC 60601-1-2: 2004/2007 (2nd/3rd Edition)	305
EMC - Electromagnetic compatibility as per IEC 60601-1-2: 2014 (4th Edition)	309
Environmental conditions	311
Regulatory information	312

Order data

Order data for country-specific cables

Description	Order number
Europe	000000-0135-487
USA	397902-9066-000
UK	309850-9003-000
Switzerland	000000-0426-193
Argentina	000000-0434-527
China	000000-0475-509
Brazil	000000-0527-729

Accessories

- Only use the accessories included in the delivery package.
- Use the accessories included in the delivery package only for the intended purpose.
- The accessories listed below are not available in all countries. Please contact your local sales agent for information on country approvals.
- Please send your orders to your local contact of the Carl Zeiss subsidiary in your country.

Designation	
INTRABEAM 600 Mammo	10MAMMOB
INTRABEAM 600 Mammo Advanced	10MAMMOA
INTRABEAM 600 Neuro	10NEUROS
INTRABEAM 600 Superficial	10SUPERF
INTRABEAM 600 GI/CR	10GICORE
INTRABEAM 600 H/N	10HEANEC
INTRABEAM 600 Upgrade	01IB600U
INTRABEAM 600 Additional Site	01IB600A
INTRABEAM NC32 Floor Stand	10FSNC32
Water phantom	60IBWAPH
Miniaturized Linear Accelerator XRS 4	60MAXRS4
Sterilization tray	20IBSTTR
Spherical applicator (Ø 1.5 cm)	20IBSP15
Spherical applicator (Ø 2.0 cm)	20IBSP20
Spherical applicator (Ø 2.5 cm)	20IBSP25
Spherical applicator (Ø 3.0 cm)	20IBSP30
Spherical applicator (Ø 3.5 cm)	20IBSP35
Spherical applicator (Ø 4.0 cm)	20IBSP40
Spherical applicator (Ø 4.5 cm)	20IBSP45
Spherical applicator (Ø 5.0 cm)	20IBSP50
Flat applicator (Ø 1.0 cm)	20IBFL10

Designation	Sales code
Flat applicator (Ø 2.0 cm)	20IBFL20
Flat applicator (Ø 3.0 cm)	20IBFL30
Flat applicator (Ø 4.0 cm)	20IBFL40
Flat applicator (Ø 5.0 cm)	20IBFL50
Flat applicator (Ø 6.0 cm)	20IBFL60
Surface applicator (Ø 1.0 cm)	20IBSU10
Surface applicator (Ø 2.0 cm)	20IBSU20
Surface applicator (Ø 3.0 cm)	20IBSU30
Surface applicator (Ø 4.0 cm)	20IBSU40
Needle applicator	20IBNE01
Sterile drapes	20IBSTDR
Radiation shields	20IBRASH
TARGIT Academy	20TARGIT

Spare parts

Using spare parts


- Use only the spare parts included in the scope of supply.
- Use the spare parts included in the scope of supply only for the described application purpose.
- Please send your orders to your local contact of the ZEISS subsidiary in your country.

Technical data

Technical data of the INTRABEAM Workplace

Mechanical variables

Operating mode	The device is designed for continuous operation.
Product lifetime	The product is designed for 7 years use
Total weight with full configuration	Max. 155 kg (341.72 lb)
Dimensions	

Electrical variables

Variable	Property
Nominal voltage ranges, switchable	(115V) 110V - 125V (230V) 220V - 240V
Rated frequency	50Hz - 60Hz
Connected load	300VA
The isolating transformer is designed for continuous operation.	

Only connect the workplace to wall outlets which are equipped with a properly connected protective ground conductor.

The workplace has a connector for potential equalization.

The line-powered components of your INTRABEAM 600 (CC 600, computer and UNIDOS E dosimeter) are run via the integrated isolating transformer (type REOMED 600).

Technical data of isolating transformer

The line-powered components in the workplace (CC 600, computer and dosimeter) are run via the integrated isolating transformer (type REOMED 600).

Variable	Property
Power connection	The isolating transformer may be connected only to wall outlets that are equipped with a properly functioning ground conductor.
Rated voltage	115 V / 230 V, switchable
Rated frequency	50 Hz - 60 Hz
Output voltage	230 VAC, 50 Hz / 60 Hz
Power consumption	140 VA, system-related
The isolating transformer is designed for continuous operation.	

Technical data of the XRS 4

Mechanics

Weight	1.62 kg
Dimensions	(175 x 110 x 70) mm
Probe	Dia. 3.2 mm, length 100 mm, coated with chromium nitride (CrN)

Radiation emission

Radiation properties

All values specified apply to operation at the maximum settings of 50 kV and 40 μ A.

10	9.8	1.5
20	2.3	6.5
30	0.9	16.7
40	0.4	37.5
50	0.25	60

Stability of radiation emission

$\pm 5\%$ for the first 5 minutes of X-ray beam emission; then $\pm 2\%$.

Radiation distribution

- Point source emission
- Isotropic distribution around isocenter of XRS
- Strong dose drop-off (approx. $1/r^3$) in soft tissue etc.
- Accuracy of the emitted dose is better than ± 1 mm with a treatment diameter of 40 mm.

Inherent filtration

Acceleration voltage [kV]	XRS HVL* [mm Al] (at water depth of 30 mm)	Effective energy [keV]
50	0.64	20.4
40	0.48	19.1

* HVL = half-value layer

Radiation exposure

The radiation exposure caused by an exposed XRS 4 at 50 kV and 40 μ A and a distance of 1 m is 9 R/h. (1.32 mGy/min).

The worst possible leakage radiation from a properly shielded XRS 4 at 50 kV and 40 μ A and at a distance of 50 mm is 0.6 mR/h. (0.0035 mGy/h).

At a distance of 50 mm, IEC 60601-2-8 permits radiation exposure amounting to ≤ 23 mR/h. (≤ 0.134 mGy/h.).

Technical data for UNIDOS E

Size	Feature	
Weight	3.2 kg (7.05 lb)	
Dimensions (W x H x D)	250 mm X 100 mm X 260 mm (9.84" X 3.94" X 10.24")	
Start-up time	5 min	
Power rating	Input voltage	(100 - 230) V $\pm$ 15%
	Rated frequency	(50 - 60) Hz
	Current consumption	0.05 A

Technical data for Control Console 600 (CC 600)

Weight	4.5 kg (9.92 lb)
Dimensions (W x H x D)	(305 x 89 x 381) mm / (12.00 x 3.50 x 15.00)"
Power connection	Via workplace isolating transformer
Input voltage	100 V - 240 V
Rated frequency	50 Hz - 60 Hz
Power consumption	60 VA
Selectable operating parameters:	
Acceleration voltage / beam current	40 kV / 40 μ A 50 kV / 5, 10, 20 or 40 μ A
	40 μ A is the maximum current of the X-ray tube over the entire voltage range.

The device is designed for continuous operation (standby).

Since the expected service life of electronic storage modules is approx. 7 years, it may be necessary to replace electronic components after this period.

Technical data of computer

Weight	12 kg (2.20 lb)
Dimensions (W x H x D)	(400 x 341 x 100) mm / (15.75 x 13.43 x 3.94)"
Power connection	Via workplace isolating transformer
Input voltage	100 V - 240 V
Rated frequency	50 Hz - 60 Hz
Power consumption	Max. 150 VA

Ethernet connection (optional),

Prerequisites for integrating the device in an existing IT network

- **When integrating the device in an existing IT network, observe the following data and measures per IEC 60601-1:**
 - The device can be connected to an existing IT network in order to transmit device parameters.
 - The IT network in which the device is to be integrated must have the following required characteristics:
 - IPv4, statically or dynamically addressed
 - The IP address space is random.
 - The IT network in which the device is to be integrated must have the following required configuration:
 - The firewall is open for the incoming ports.
 - Transmissions from the device to the outside are not planned.
 - The network connection with the device is prescribed as follows:
 - The IT network has an Ethernet connection via RJ45, at least Cat. 5, 100 Mbit/s.
 - The device has a network isolator as per IEC 60601-1, and is therefore electrically separated from the network.
 - Information flow and routing between the device, the IT network and other devices:
 - The essential performance features of the device do not require any connection to the IT network or data transmission.
 - The topology of the IT network and the required routing lie are the responsibility of the operator.

If the IT network cannot provide the features required to serve the purpose of integrating the device in the IT network, the following hazardous situations may result:

- Increased earth leakage current to the power output sockets:
Danger of an electrical shock if the user contacts a live housing part and touches the patient.

- Virus infestation of the system:
Risk for the patient due to overexposure and/or longer operation and anesthesia times.

Integration of the device in an IT network that includes other ME devices can lead to risks which were not previously known to exist for patients, operators or third persons.

- The operator of the device must determine, analyze, evaluate and manage these risks before integrating the device in the IT network.
- The following changes in the IT network can lead to new risks, and therefore require additional analyses:
 - Changes in the IT network configuration
 - Connecting additional elements to the IT network
 - Removing elements from the IT network
 - Update of devices connected to the IT network
 - Upgrade of devices connected to the IT network

Data carriers must be virus-free.

The operator is responsible for network security and ensuring that the data carriers used (USB sticks, CDs and DVDs) are free of viruses.

Virus-infested software creates risks for the patient due to overexposure and/or longer operation and anesthesia times.

External interlock and X-radiation indicator

WARNING


Danger due to radiation during treatment

During irradiation, X-ray radiation can emerge from the patient into his surroundings, thus harming the user and other persons in the treatment room.

- It is essential for you to equip your treatment room with an external interlock (interlocking device) and an X-ray display (warning lamp).
- Before beginning irradiation, ensure that all persons who are not wearing suitable radiation protection suits have left the treatment room.

The external interlock connector "EXT INTERLOCK" located on the back of the workplace has the following two functions:

- Option of connecting a separate switch to stop, prevent or release the generation of X-ray radiation.
- Option of connecting an additional external display (warning lamp) to indicate that X-ray radiation is being generated.

The interlock plug connected to the workplace fits the "EXT INTERLOCK" connector. Pins 1 and 2 of this interlock connector are connected by a wire jumper. If this interlock plug is plugged into the "EXT INTERLOCK" connector, it completes the external interlock circuit.

Adjacent illustration: Pin assignment of "EXT INTERLOCK" connector.

If you would like to use a separate switch, for example a door switch, you will have to connect it to pin 1 and 2 on the "EXT INTERLOCK" connector using a separate plug. This connection will be implemented by the customer.

A suitable plug for the external interlock connector is supplied with the control unit. You can also purchase the appropriate plug from ZEISS Service as a spare part, see page 292 Spare Parts "Interlock Control Plug".

NOTE**Wrong or missing connection to external interlock**


If mains power is supplied to the workplace by mistake via the external interlock, this will cause a fault in the INTRABEAM 600. If the user has not established a proper cable connection to the external interlock, no treatment is possible.

- Make sure that the external interlock is properly connected to the workplace.

Connecting the external switching contact for the interlock

Use a low-impedance connection (< 10 ohms) between pin 1 and pin 2 of the plug connector.

Circuit recommended for connection to pins 1 and 2:

External X-ray display

The external X-ray display is supplied by the customer.


ZEISS recommends using an ohmic resistance (R) on pins 3 and 4 for current limitation and a relay for controlling a lamp of the external X-ray display.

If X-ray radiation is generated, pins 3 and 4 form a closed circuit; otherwise they form an open circuit.

Maximum voltage on pins 3 and 4 50 V DC or 24 V AC

Maximum current on pins 3 and 4 0.5 A rms

Example:

The connection possibility explained above is only a basic principle. Coupling is possible only if the requirements from IEC 60601-1-1: 2000 or IEC 60601-1:2005, Chapter 16 (Safety of medical electrical systems) are met and a proper installation is performed. Furthermore, it must be ensured that the voltage source and further separations between the line voltage and low voltage are carried out according to IEC 60601-1. For the installation, the pertaining national regulations must be observed.

EMC (electromagnetic compatibility)

Medical electrical equipment and medical systems are subject to special precautions with regard to electromagnetic compatibility (EMC).

In order to avoid the occurrence of EMC interference, the INTRABEAM 600 may only be installed, started up and maintained in the manner indicated in these Instructions for Use and only with components supplied by ZEISS.

- Please follow the EMC guidelines on the following pages.

CAUTION

Danger from electromagnetic radiation

Electrical devices can influence each other due to their electromagnetic radiation. The use of non-approved components (accessories, transformers of all types, cables) can cause increased emissions or reduce the device's immunity.

- The INTRABEAM 600 may be installed and operated in the immediate vicinity of another electrical device only in the device combinations described in these Instructions for Use.
If it is unavoidable to operate the INTRABEAM 600 in the vicinity of other electrical devices, proper functioning of the INTRABEAM 600 must be monitored.
- Only use accessories, transformers, cables and spare parts which are specified in these Instructions for Use or which have been approved by ZEISS for the INTRABEAM 600.
- Do not use any portable or mobile RF communication equipment or transmitters, including peripheral devices such as antenna cables or external antennas, in the proximity of the INTRABEAM 600 (minimum distance: 30 cm).
It cannot be ruled out that the functioning of the INTRABEAM 600 will be influenced or its performance will be impaired as a result of this.

Medical electrical systems

Any additional equipment connected to medical electrical devices must demonstrably comply with the applicable standards and directives (e.g. IEC 60950-1:2005, modified + Cor.:2006 + A1:2009, modified + A1:2009/Cor.:2012 + A2:2013, modified for data processing equipment). Furthermore, any configurations must comply with normative requirements for electronic medical devices (see IEC 60601-1:2005 + Cor.:2006 + Cor.:2007 + A1:2012).

Anyone connecting additional equipment to medical electrical devices is a system configurer and as such is responsible for compliance of the system with the normative requirements for medical electrical systems.

Please note that local legislation takes priority over the above-mentioned normative requirements. If you have any questions, please contact your local dealer or ZEISS Service.

EMC-relevant full configuration of the INTRABEAM 600

The electromagnetic compatibility was determined on a full configuration of the INTRABEAM 600 with XRS 4.

The full configuration of an INTRABEAM 600 with XRS 4 comprises the following components:

- INTRABEAM Workplace
(Standard package, sales code: 01PIBCONF)
- Stand (NC32 INTRABEAM Floor Stand)
Order number: 000000-1096-945
- XRS 4 on stand

The full configuration contains the following exposed cable connections:

- Power supply for workplace: power cord, length 6 m
- Control cable for XRS 4, length 15 m
- Control cable for PDA/PAICH, length 1.5 m
- Control cable for ionization chamber ("Ion Chamber" connector)
- Power supply for stand: power cord, length 6 m

WARNING

If you use any components other than the ones approved by ZEISS for the INTRABEAM 600 (i.e. ME devices, accessories, isolating transformers (converters) and cables) without express permission, this may lead to increased electromagnetic radiation or a reduced electromagnetic immunity of your INTRABEAM 600.

Safe, trouble-free operation of your INTRABEAM 600 is then no longer possible, which may lead to overexposure of the patient.

- Only use the components approved by ZEISS in your INTRABEAM 600.

EMC - Electromagnetic compatibility as per IEC 60601-1-2: 2004/2007 (2nd/3rd Edition)

Electromagnetic interference emissions

Table 1: The INTRABEAM 600 is intended for operation in an electromagnetic environment as specified below. The operator and the user must ensure that the INTRABEAM 600 is used only in such an environment.

RF emissions as per CISPR 11	Group 1	The INTRABEAM 600 uses RF energy only for its internal functions. As a result, RF emissions are very low and unlikely to cause any interference in nearby electronic devices.
RF emissions as per CISPR 11	Class B	The INTRABEAM 600 is suitable for use in all facilities including locations in residential environments and those directly connected to the public power supply network which also supplies buildings used for residential purposes.
Harmonic emissions as per EN 61000-3-2	Class A	
Emissions of voltage fluctuations/ flicker as per EN 61000-3-3	Compliant	

Electromagnetic immunity

Table 2: The INTRABEAM 600 is intended for operation in an electromagnetic environment as specified below. The operator and the user must ensure that the INTRABEAM 600 is used only in such an environment.

Electrostatic discharge (ESD) as per IEC 61000-4-2	±6 kV contact discharge	±6 kV contact discharge	Floors should be made of wood or cement or covered with ceramic tiles. If the flooring contains synthetic materials, the relative humidity must be at least 30%.
	±8 kV air discharge	±8 kV air discharge	

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment guidelines
Fast transient/burst immunity as per IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input and output lines	±2 kV for power supply lines ±1 kV for input and output lines	The quality of the supply voltage should be that of a typical business or hospital environment.
Surges as per IEC 61000-4-5	± 1 kV phase-neutral voltage ± 2 kV phase/neutral to ground voltage	± 1 kV phase-neutral voltage ± 2 kV phase/neutral to ground voltage	The quality of the supply voltage should be that of a typical business or hospital environment.
Voltage dips, short interruptions and supply voltage fluctuations as per IEC 61000-4-11	< 5 % U_T (> 95 % dip in U_T) for 1/2 cycle 40 % U_T (60 % dip of U_T) for 5 cycles 70 % U_T (30 % dip of U_T) for 25 cycles < 5 % U_T (95 % dip in U_T) for 5 s	< 5 % U_T (> 95 % dip in U_T) for 1/2 cycle 40 % U_T (60 % dip of U_T) for 5 cycles 70 % U_T (30 % dip of U_T) for 25 cycles < 5 % U_T (95 % dip in U_T) for 5 s	The quality of the supply voltage should be that of a typical business or hospital environment. If the user of the INTRABEAM 600 requires continued operation even in the event of interruptions in the power supply, we recommend powering the INTRABEAM 600 from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field as per IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields in the supply frequency should correspond to the typical values that are found in business and hospital environments.

NOTE: U_T is the AC supply voltage prior to application of the test level.

Electromagnetic immunity for ME equipment or systems that are not life-sustaining

Table 3 The INTRABEAM 600 is intended for operation in an electromagnetic environment as specified below. The operator and the user must ensure that the INTRABEAM 600 is used only in such an environment.

Portable and mobile radio communication equipment should not be used closer to the INTRABEAM 600, including cables, than the recommended safety distance that is calculated using the equation applicable to the transmission frequency involved.

Conducted RF disturbances as per EN 61000-4-6	3 V _{effective value} 150 kHz to 80 MHz	3V	Recommended separation distance:	
			$d = 12\sqrt{P}$	
Radiated RF disturbances as per EN 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 12\sqrt{P}$	for 80 MHz to 800 MHz
			$d = 23\sqrt{P}$	for 800 MHz to 2.5 GHz

where P is the output power rating of the transmitter in watts (W) according to the transmitter manufacturer's specifications and d is the recommended safety distance in meters (m).

The field strength of stationary transmitters as determined by a site survey should be^a less than the compliance level in all frequency ranges.^b

Interference may occur in the vicinity of equipment marked with the following symbol.

Note 1 At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is influenced by absorption and reflection by structures, objects and persons.

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment Guidelines
----------------	----------------------	------------------	--

^a Theoretically, field strengths of stationary transmitters such as base stations for mobile telephones and mobile land radio equipment, amateur radio stations, AM and FM radio broadcast and TV broadcast transmitters cannot be predicted accurately. To assess the electromagnetic environment with respect to stationary RF transmitters, a site study of the electromagnetic phenomena should be considered. If the field strength measured at the location where the device is used exceeds the above compliance levels, the device must be monitored to verify proper functioning. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

^b Over the frequency range of 150 kHz to 80 MHz, the field strength should be less than 3 V/m.

Recommended safety distances

Table 4 The INTRABEAM 600 is intended for use in an electromagnetic environment in which RF disturbances are controlled. The operator and the user of the INTRABEAM 600 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the INTRABEAM 600 - depending on the output power of the communication equipment as specified below.

Rated power of transmitter [W]	Separation distance depending on transmission frequency [m]		
	150 KHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 12\sqrt{P}$	$d = 12\sqrt{P}$	$d = 23\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation indicated for each column, with P being the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer's specifications.

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is influenced by absorption and reflection by structures, objects and persons.

EMC - Electromagnetic compatibility as per IEC 60601-1-2: 2014 (4th Edition)

Electromagnetic interference emissions

The INTRABEAM 600 is intended for operation in an electromagnetic environment as specified below. The operator and the user must ensure that the INTRABEAM 600 is used only in such an environment.

RF emissions according to CISPR 11	Group 1
RF emissions according to CISPR 11	Class B
Emissions of harmonics as per IEC 61000-3-2	Class A
Emissions of voltage fluctuations/flicker as per IEC 61000-3-3	Compliant

Electromagnetic immunity

The INTRABEAM 600 is intended to be operated in the electromagnetic environment specified below. The operator and the user must ensure that the INTRABEAM 600 is used only in such an environment.

Electrostatic discharge (ESD) as per IEC 61000-4-2	±8 kV contact discharge	±8 kV contact discharge
	±15 kV air discharge	±15 kV air discharge
Fast transient/burst immunity as per IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines
	±1 kV for input and output lines	±1 kV for input and output lines
Surges as per IEC 61000-4-5	±1 kV phase-neutral voltage	±1 kV phase-neutral voltage
	±2 kV phase/neutral to ground voltage	±2 kV phase/neutral to ground voltage

Immunity tests	IEC 60601 - test level	Compliance level
Magnetic field for supply frequency (50/60 Hz) as per IEC 61000-4-8	30 A/m	30 A/m
Voltage drops, short interruptions and supply voltage fluctuations as per IEC 61000-4-11	0 % U_T for 1/2 cycle	0 % U_T for 1/2 cycle
	0 % U_T for 1 cycle	0 % U_T for 1 cycle
	70 % U_T for 25/30 cycles	70 % U_T for 25/30 cycles
	0 % U_T for 250/300 cycles	0 % U_T for 250/300 cycles

Electromagnetic immunity for ME equipment or systems that are not life-sustaining

The INTRABEAM 600 is intended to be operated in one of the electromagnetic environments specified below. The operator and the user must ensure that the INTRABEAM 600 is used only in such an environment.

Immunity tests	IEC 60601 - test level	Compliance level
Conducted RF disturbances as per EN 61000-4-6	3 V 150 kHz to 80 MHz	3V
	6 V ISM bands between 150 kHz and 80 MHz	6V
Radiated RF disturbances as per EN 61000-4-3	3 V/m 80 MHz - 2.7 GHz	3 V/m
Radiated RF disturbances by electromagnetic near fields of wireless communication systems as per EN 61000-4-3	27 V/m 385 MHz	27 V/m
	28 V/m 450 MHz, 810 MHz - 2.45 GHz	28 V/m
	9 V/m 710 MHz - 780 MHz, 5.24 GHz - 5.785 GHz	9 V/m

Environmental conditions

For operation	Temperature	+ 15°C ... + 38°C + 59°F ... + 100°F
	Relative humidity	30% ... 75%
	Atmospheric pressure	800 hPa ... 1060 hPa
For transport and storage	Temperature	- 20 °C ... + 60°C - 4°F ... + 140°F
	Relative humidity (without condensation)	10%...85%
	Atmospheric pressure	600 hPa...1060 hPa

If the UNIDOS-E is shipped separately from the complete INTRABEAM 600 device, the following conditions shall apply:

UNIDOS E, for transport and storage	Temperature	- 20 °C ... + 60°C - 4°F ... 140°F
	Relative humidity (without condensation)	10% ... 85%
	Atmospheric pressure	600 hPa ... 1200 hPa

Regulatory information

Description	Labeling
Electrical safety	The medical device complies with the requirements according to – IEC 60601-1:2005 + Cor. :2006 + Cor. :2007 + A1:2012, – IEC 601-1:1988+A1:1991+A2:1995, – CAN/CSA C22.2 60601-1:2014. Classification according to the level of protection against electric shock: Protection class 1 Protection type (according to IEC 60529): PI X0
Electromagnetic compatibility	The medical device is classified according to IEC 60601-1-2 of Group 1, Class B (as per CISPR 11).
Classification	The medical device is classified according to Annex IX of Medical Device Directive 93/42/EEC as Class IIb.
CE marking	The medical device meets the essential requirements stipulated in Annex I of Medical Device Directive 93/42/EEC. In line with Annex XII of Medical Device Directive 93/42/EEC, the medical device is labeled with 0297.
CSA mark	The medical device is marked with .

Supplementary documents

Applicable user manuals314

Literature on the INTRABEAM 600315

Checklist for OR equipment316

Applicable user manuals

- NC32 INTRABEAM Floor Stand
ZEISS order number: 1096-945
Instructions for Use, ZEISS: G-30-1476
- INTRABEAM Drape
ZEISS order number: 32090-0000-000
Instructions for Use, ZEISS: G-30-1526
- Radiation shield, flat
ZEISS order number: 304534-6000-595
Instructions for Use, ZEISS: G-30-1890
- INTRABEAM Spherical Applicator 1.5 cm, 2.0 cm, 2.5 cm, 3.0 cm
3.5 cm, 4.0 cm, 4.5 cm, 5.0 cm and Set
ZEISS order numbers: 304534-6000-557, -558, -559, -560, -561,
-562, -563, -564 and -348 (set)
Instructions for Use, ZEISS: G-30-1478
- INTRABEAM Needle Applicator
ZEISS order number: 304534-3001-200
Instructions for Use, ZEISS: G-30-1723
- INTRABEAM Surface Applicator Set 1.0 cm, 2.0 cm, 3.0 cm and 4.0 cm
ZEISS order numbers: 304534-3001-530, -531, -532, -533
Instructions for Use, ZEISS: G-30-1858
- INTRABEAM Flat Applicator Set 1.0 cm, 2.0 cm, 3.0 cm, 4.0 cm,
5.0 cm and 6.0 cm
ZEISS order numbers: 304534-3001-565, -566, -567, -568, -569 and -
570
Instructions for Use, ZEISS: G-30-1880
- UNIDOS E dosimeter from PTW, Type 10010
PTW Instructions for Use, English: D545.131.00/06 en
PTW Instructions for Use, German: D545.151.00/06 de
- INTRABEAM Water Phantom
ZEISS order number: 304534-7000-600
Instructions for Use, ZEISS: G-30-1767.

Literature on the INTRABEAM 600

Current clinical, physical and scientific publications on the INTRABEAM 600 and the respective fields of application can be found under www.zeiss.com/radiotherapy

Checklist for OR equipment

Items of equipment according to specification:	Check
Workplace with power cable	
Workplace equipped with CC 600, computer, keyboard, touch-screen and dosimeter	
Equipment for the X-ray beam source: (transport and storage in the workplace)	
– XRS 4 with probe protector	
– XRS cable	
– (CC 600 integrated in the workplace)	
Accessories for X-ray beam source: (transport and storage in the workplace)	
– PAICH (probe adjuster and ionization chamber holder)	
– PDA (photo diode array)	
– PAICH/PDA cable	
Equipment for ionization chamber (IC): (dosimeter in the workplace)	
– IC (ionization chamber)	
– IC extension cable	
– IC shield	
Sterile applicators	
Instructions for Use	
Radiation monitoring device (not a component of the INTRABEAM 600 - recommendation)	
Scale, sterile (not a component of the INTRABEAM 600 - recommendation)	

Index

A

About menu	138
Acceleration voltage, CC 600	297
Accessories	291
Accessories for inspection	16, 64
Accessories for the INTRABEAM 600	18, 65
Accessories for your INTRABEAM 600	86
Address of manufacturer	2
Administering XRS X-ray sources	142, 152
Administrative tools	125
Administrative tools, user administration	115
Administrator	115
Air pressure sensor	69
Air pressure sensor, calibration	264
Anesthetics	43
Applicators	88
Applicators, Equipment inventory	152
Archiving, non-permanently	36
Archiving, not permanent	36
Arrangement of equipment and personnel	99

B

Beam current, CC 600	297
Beam deflection adjustment	173
Beam deflection, window	218
Before every use	32
Beryllium (Be)	74
Breast cancer	100

C

Calibrating the IC ionization chamber	255
Calibrating the temperature sensor	267
Calibrating the XRS 4	254
Calibration	257
Calibration instructions	254
Calibration, screen mask	257
Calibrations	39
Calibrations at the manufacturer's plant	257
Canceling a treatment	224

CC 600 Control Console	64, 111
CC 600, functionality of	69
Centering (straightening), XRS probe	80
Changing the line fuse	277
Checking line connections of devices	97
Checklist	210
Checklist for OR equipment	316
Class B (practice area)	28
Cleaning	44, 279, 281
Cleaning and disinfecting components	282
Cleaning components	281
Cleaning the touchscreen	283
Clinical consideration	39, 40, 213, 219
Components for inspection	16, 64
Components for quality assurance	64
Components, approved	97
Components, line-powered	97
Computer of the INTRABEAM 600	72
Computer shutdown	33
Computer software login	112
Computer, technical data	297
Conducting treatment	212
Connecting	102
Continue with previous patient	225
Continuous operation	294
Continuous operation (standby)	69
Continuous operation (standby), CC 600	297
Control Console (CC 600)	69
Control console (CC 600)	14
Control Console, CC 600	111
Control lamps	38, 69
Control panel display	70
Control panel of the CC 600	70
Controls, displays, connections	70
Conversion of the system	32
Counting rate of external X-ray monitor	218
Counting rate of internal X-ray monitor	218
Cross hair	200
Current consumption	296

D

Damage and modifications, obvious	97
Data carriers must be free from viruses.	36
Data loss	36
Data Management	120
Data protection	35

Depth dose curve	200
Depth dose curve diagram	200
Depth dose curve, measurement of	89
Designation in these Instructions for Use	3
Determining the beam centering	258
Determining the linearity of the IRM internal X-ray monitor	262
Dimensions	297
Dimensions, CC 600	297
Dimensions, equipment cart	293
Dimensions, UNIDOS E	296
Disinfectants	284
Disinfectants, commercial	284
Disinfecting	44, 279, 281
Disinfecting surfaces	284
Disk space, no longer sufficient	36
Display for setting parameters	69
Displaying the treatment record	242
Disposal	288
Disposal of components of the INTRABEAM 600	45
Disposal of consumables	45
Disposal of disposable components	45
Dosimeter UNIDOS E	84
Drape	3
Dynamic Offsets	173

E

Editing folders	122
Electrical requirement	94
Electrical safety	31, 34
Electrics of the INTRABEAM 600	40
Electromagnetic compatibility	28, 303
Electromagnetic interference	42
Electromagnetic radiation	28
EMC - Electromagnetic compatibility	305
EMC - electromagnetic compatibility	309
EMC (electromagnetic compatibility)	303
EMC requirements of IEC 60601-1-2	28, 42
EMC, respective precautions	42
Emergency backup power supply	31
Emergency Off button	70, 224
EMERGENCY STOP	38
Environmental conditions	28
Equipment inventory	140
Equipment inventory Applicators group	152
Equipment inventory, "XRS linear accelerator" group	142
Equipment inventory, Equipment group	148

Equipment preparation	278
Equipment, Equipment inventory	148
Ethernet connection (optional)	298
Exporting files	122
External interlock	300

F

Field of application	24
File, menu	124
Fire	43
Flowchart	105, 159
Fractional treatment, generating a treatment record	244

G

General plan	99
Generating a short form treatment record	248
Generating a treatment record	244
Generating a treatment record of a fractional treatment	244, 245
Generating a treatment record of a single treatment	250
Generating X-ray beams	37
GMV company, Spain	2

H

Handling	278
Handling and preparing the medical device	24
Handling components	278
Handling, transport and shipment	278
Hazard symbols	23
High voltage generator	74

I

Importing calibration data	119
Indications for use	24
Information for the operator	27
Information security	35
Information symbols	23
Information, required additional	30
Initial calibrations	257
Input voltage	297
Input voltage, CC 600	297
Input voltage, UNIDOS E	296
Inspection instruments	78
Instructions for Use, present	30
Intended use	24

Interference with device, unauthorized	29
Interlock connector, external	300
Interlock status	37
Interlock system, optical	37
Interlock, external	38
Interrupting the treatment	222
INTRABEAM 600	14
INTRABEAM 600 software	14
INTRABEAM 600 user interface	14
INTRABEAM 600, checklist for the	210
INTRABEAM 600, disposal of components	45
INTRABEAM 600, electrics of the	40
INTRABEAM 600, safety devices	37
INTRABEAM 600, use of	34
INTRABEAM Flat Applicator Set	30, 314
INTRABEAM Needle Applicator	30, 314
INTRABEAM Spherical Applicators	30, 314
INTRABEAM Start Manager	113
INTRABEAM Surface Applicator Set	30, 314
INTRABEAM Water Phantom	30, 314
INTRABEAM Workplace	66
INTRABEAM Workplace, mobile workstation	64
Ionization chamber holder, PAICH	80
IRM internal X-ray monitor	64
IRM linearity	262
IRM X-ray monitor, determining linearity	262
IRM X-ray monitor, internal	74, 75
Irradiation	39
Irradiation dose, delivered	226
Irradiation dose, specified	226
Isolating transformer	97
Isolating transformer in equipment cart	294
Isotropy	182
Isotropy check	176
Isotropy of probe tip	83
Isotropy, adjustment of	182
Isotropy, measurement of	89

K

Key to symbols	23
----------------------	----

L

Labeling of the XRS 4	56
Legal regulations	27
Liability and warranty	29

Limit class A	42
LINAC Configuration tool	118
Line voltage, required	94
Line-powered components	294
Loss of warranty	29

M

Main monitoring system	75
Major performance characteristics	25
Malfunction	34
Manufacturer's address	2
Marking the treatment site	192
Medical electrical systems	27, 303
Memory space, insufficient	36
Menu About	138
Menu bar	108
Menu, File	124
Menu, QA Reports	133
Menu, Tools	125
Modifications, do not make any	29
Multiple sockets, none	97

N

NC32 INTRABEAM Floor Stand	30, 314
Needle applicator	44, 280
Nominal voltage range, equipment cart	294
Nonsterile area	99
Normal use	25

O

Obligations of the operator	27
"Observe transport position" warning label	50, 95
Offline mode	116
Online mode	116
Operating mode, equipment cart	293
Operation	209
Operation of the device, termination	33
Optical interlock system	34, 37
OR equipment, checklist for	316
Order data	290
Order data for country-specific cables	290
Output check of dose rate.	182
Output voltage, isolating transformer in equipment cart	294
Oxygen	43

P

PAICH	80
Password	115
Password entry	205
Pause button	218
Pause mode	222
PDA XRS check	182
PDA, photo diode array	83
Performing quality assurance	160
Photo diode array (PDA)	83
Physician, radiation therapist	205
Physicist	205
Planning fractional treatment	202
Potential equalization	38
Power connection	297
Power connection, CC 600	297
Power consumption	297
Power consumption, CC 600	297
Power consumption, maximum	94
Power switch	111
Precautions regarding EMC	42
Preparations for use	91
Preparing and handling the medical device	24
Preparing the workplace for use	94
Principal components	69
Privacy protection	35
Probe adjuster	80
Probe protector	77
Probe tip	74
Product information	30
Product lifetime, equipment cart	293
Proper handling of components during quality assurance	164
Protection components	77
Protection of components of INTRABEAM 600	64
Protective clothing	34
Protective ground contact	32
Protective ground contact, correctly connected	94
PTW company	2

Q

QA Reports menu	133
Qualification of the user, further	24
Quality assurance	158
Quality assurance with PDA	83

R

radiance software	117, 118
Radiation distribution	295
Radiation emission, stability of	295
Radiation exposure	296
Radiation output, monitored	39
Radiation pattern, aspherical	75
Radiation properties	295
Radiation shielding	34
Radiation therapist, treating physician	205
Radio frequency generating devices	28, 42
Rated frequency	297
Rated frequency, CC 600	297
Rated frequency, UNIDOS E	296
Records, treatment and verification	69
Reference to calibration	34
Registered trademark of the GMV company	2
Registered trademark of the PTW company	2
Regulatory information	312
Reports menu	181, 188
Requirements for system operation	93
Responsibilities for data protection	35
Review Treatment	231
"Risk of crushing" warning label	50

S

Safety devices of the INTRABEAM 600	37
Safety during tests and verification	39
Safety during treatment	33, 39
Safety measures	21
Saving the calibration file	146
Select working status	116
Selecting a patient	190
Selecting an XRS X-ray source	196
Service Identification Program (SIP)	271
Service requirements	271
Shipment	278
Signing the treatment record personally	205
SIP label	50
SIP label with SIP number	271
Software and electronics	69
Software user interface	106
Solid water, water-equivalent material	81
Space requirement	93
Spare parts	292

Specialists, qualified	24
Specialists, trained and briefed	24
Start button	218
Start of INTRABEAM software	114
Starting the computer	112
Starting the INTRABEAM software	114
Starting the software	111
Starting the treatment	218
Startup of the INTRABEAM 600	111
Start-up time, UNIDOS E	296
Sterile area	99
Sterilization	44, 279
Sterilization for interstitial treatment	44, 279
Sterilization for tumor bed treatment	44, 279
Sterilizing	44, 279
STOP button	224
Straightness measurement, XRS probe	80
Structure of the software user interface	106
Supplementary documents	313
Surface temperature of applicators	40, 213, 219
Switching the workplace on	111
System check, printing reports	181, 188
System description	64

T

Tab bar	107
Technical data	293
Technical data for Control Console 600	297
Technical data of computer	297
Technical data, UNIDOS E dosimeter	296
Technical data, XRS 4	295
Tests and verification, safety during	39
Thermal energy, generation of	40, 213, 219
Tip adjustment with PAICH	168
Tools, menu	125
Trademark	2
Transport	31, 278
Treatment approval	205
Treatment approval screen mask	205
Treatment approval, screen mask	205
Treatment facility	99
Treatment facility, spatial arrangement	99
Treatment parameters, display	231
Treatment prescription	190
Treatment record, generating short form	248
Treatment room	99

Treatment terminated	226
Treatment, safety during	33, 39
Troubleshooting table	271
Troubleshooting	227
Tumor bed treatment, arrangement for	100

U

UNIDOS E dosimeter	84
UNIDOS E dosimeter from PTW	30, 314
UNIDOS E dosimeter, technical data	296
UNIDOS E, start-up time	296
Use of the INTRABEAM 600	34
User	115
User groups	115
User interface	106
User manuals, other applicable	314
User name and password	115

V

Verification devices	80
Verification equipment, administer all	148
Verification, validity of	160
V-guide	78
Video instruction, video window	166
Visual system check	97

W

Warranty and liability	29
Water phantom	89
Water-equivalent material, solid water	81
Weight	297
Weight, CC 600	297
Weight, equipment cart	293
Weight, UNIDOS E	296
Workflow	104
Workplace	66

X

X-ray beam interlock, external	301
X-ray display, external	302
X-ray exposure time	34
X-ray monitor, internal	39
X-ray probe	74
X-ray radiation, defined, isotropic, constant	25

X-ray radiation, external display	300
X-ray source, status of	142, 148
X-ray sources, equipment inventory	142
XRS 4	74
XRS 4 X-ray source	64, 74
XRS 4, labeling of	56
XRS 4, technical data	295
XRS mount, area	78
XRS tip adjustment	168
XRS X-ray source, annual calibration	146
XRS X-ray source, check	89
XRS X-ray source, symbols on	56

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page. There is no handwriting or other markings on the paper.

[illegible]

Carl Zeiss Meditec AG
Goeschwitzer Strasse 51-52
07745 Jena
Germany

Fax: + 49 (0)7364 - 20 4823
Email: info.meditec@zeiss.com
Internet: www.zeiss.com/med

Annex I

Information of the manufacturer and developer of the medical device

Developer	
Name of the legal entity with transliteration	Carl Zeiss Meditec AG
Developer location address	Rudolf-Eber-Strasse 11, 73447 Oberkochen, Germany
Manufacturer	
Name of the legal entity with transliteration	Carl Zeiss Meditec AG
Legal entity location address	Goeschwitzer Strasse 51-52, 07745 Jena, Germany
Telephone numbers	Tel.: +49 36 41/220 439, Fax: +49 36 41/220 756
Legal entity e-mail	info@meditec.zeiss.com
Medical device manufacture place	Carl Zeiss Meditec AG, Rudolf-Eber-Strasse 11, 73447 Oberkochen, Germany

Information about the authorized representative of the manufacturer in the territory of the Russian Federation

Legal organization form and full name of the legal entity	Carl Zeiss Limited Liability Company
Short name of the legal entity (when available)	Carl Zeiss LLC
Legal entity location address	Room 21, office 1, floor 4, 29 Serebryanicheskaya embankment, Moscow, 109028, Russia
Telephone numbers	Telephone (495) 933-51-51; Fax (495) 933-51-55
Legal entity e-mail	office@ru.zeiss.com

Annex 2. Definitions

Application for registration of medical device	Operational documentation, extract from the manufacturer's technical documentation
I. INTRABEAM Intraoperative Radiotherapy System, version INTRABEAM PRS 500 Consisting of:	INTRABEAM system, INTRABEAM PRS 500, system
1. X-ray source (XRS 4).	XRS source, X-ray source XRS, radiation source XRS 4, XRS 4
2. V-Block.	"V-block" guide
3. X-Block.	X-block
4. Control console PRS 500.	Control console PRS 500
5. Terminal PRS 500 (Computer with touchscreen and INTRABEAM Software).	User terminal PRS 500, terminal PRS 500
6. Probe Adjuster Ion Chamber Holder (PAICH).	Probe centering device / holder PAICH, PAICH (probe centering device and ionization chamber holder)
7. Ionization chamber / Ionization Chamber Shield.	Ionization chamber IC, Ionization chamber "IC shield"
8. Photo diode array (PDA).	Photodiode array PDA.
9. Electrometer (UNIDOS E) manufactured by PTW Freiburg Physikalisch-Technische Werkstätten Dr. Pöchlau GmbH, production site located at: Lörracher Str. 7, 79115 FREIBURG, GERMANY.	Dosimeter UNIDOS E, Electrometer UNIDOS E, UNIDOS E
10. INTRABEAM Water Phantom.	Water Phantom, INTRABEAM Water Phantom
11. Isolating Transformer.	---
12. INTRABEAM Cart.	INTRABEAM Workplace, Cart INTRABEAM
13. NC32 INTRABEAM Floor Stand.	Floor Stand NC32 INTRABEAM, Floor Stand
14. Applicators for INTRABEAM Intraoperative Radiotherapy System in sets and on a standalone basis, with accessories P3H 2014/2108 of 01.11. 2018 г. - INTRABEAM Spherical Applicator 1.5 cm (set of INTRABEAM Spherical Applicator 1.5 cm, calibration disc) - INTRABEAM Spherical Applicator 2.0 cm (set of INTRABEAM Spherical Applicator 2.0 cm, calibration disc)	INTRABEAM Spherical Applicator, Spherical Applicator

- INTRABEAM Spherical Applicator 2.5 cm (set of INTRABEAM Spherical Applicator 2.5 cm, calibration disc) - INTRABEAM Spherical Applicator 3.0 cm (set of INTRABEAM Spherical Applicator 3.0 cm, calibration disc) - INTRABEAM Spherical Applicator 3.5 cm (set of INTRABEAM Spherical Applicator 3.5 cm, calibration disc) - INTRABEAM Spherical Applicator 4.0 cm (set of INTRABEAM Spherical Applicator 4.0 cm, calibration disc) - INTRABEAM Spherical Applicator 4.5 cm (set of INTRABEAM Spherical Applicator 4.5 cm, calibration disc) - INTRABEAM Spherical Applicator 5.0 cm (set of INTRABEAM Spherical Applicator 5.0 cm, calibration disc) - INTRABEAM Flat Applicator 1.0 cm (set of INTRABEAM Flat Applicator 1.0 cm, INTRABEAM Position Marker 1.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Flat Applicator 2.0 cm (set of INTRABEAM Flat Applicator 2.0 cm, INTRABEAM Position Marker 2.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Flat Applicator 3.0 cm (set of INTRABEAM Flat Applicator 3.0 cm, INTRABEAM Position Marker 3.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Flat Applicator 4.0 cm (set of INTRABEAM Flat Applicator 4.0 cm, INTRABEAM Position Marker 4.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Flat Applicator 5.0 cm (set of INTRABEAM Flat Applicator 5.0 cm, INTRABEAM Position Marker 5.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Flat Applicator 6.0 cm (set of	INTRABEAM Flat Applicator, Flat Applicator
- INTRABEAM Flat Applicator 6.0 cm, INTRABEAM Position Marker 6.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Surface Applicator 1.0 cm (set of INTRABEAM Surface Applicator 1.0 cm, INTRABEAM Position Marker 1.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Surface Applicator 2.0 cm (set of INTRABEAM Surface Applicator 2.0 cm, INTRABEAM Position Marker 2.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Surface Applicator 3.0 cm (set of INTRABEAM Surface Applicator 3.0 cm,	INTRABEAM Surface Applicator, Surface Applicator

INTRABEAM Position Marker 3.0cm, INTRABEAM Lumen Plug, calibration disc) INTRABEAM Surface Applicator 4.0 cm (set of INTRABEAM Surface Applicator 4.0 cm, INTRABEAM Position Marker 4.0cm, INTRABEAM Lumen Plug, calibration disc)	
15. Container for sterilizing INTRABEAM Spherical Applicators (including re-usable filter for container)	Sterile container
16. Transport case	-
17. Dustcover	-
18. INTRABEAM DRAPE, place of production Microtek Dominicana, S.A., production site located at: Zona Franca II, La Romana 22000, Dominican Republic	-
19. Portable data storage device (USB 3.0 Memory stick, 32GB)	-
20. XRS Cable	XRS cable
21. Network/mains cable	
22. PAICH/PDA Cable	-
23. Operational documentation including: 23.1 INTRABEAM PRS 500 Instruction for use 23.2 INTRABEAM Spherical Applicator Instruction for use 23.3 INTRABEAM Surface Applicator Instruction for use 23.4 INTRABEAM Flat Applicator Instruction for use 23.5 INTRABEAM Water Phantom Instruction for use 23.6 NC32 INTRABEAM Floor Stand Instruction for use 23.7 UNIDOS E Instruction for use	
II. INTRABEAM Intraoperative Radiotherapy System, version INTRABEAM 600 Consisting of:	INTRABEAM system, INTRABEAM 600, system
1. X-ray source (XRS 4).	XRS source, X-ray source XRS, radiation source XRS 4, XRS 4
2. Touchscreen	Monitor (with touch screen)
3. Control console 600 (CC600)	Control console CC 600 (Control unit CC 600)
4. INTRABEAM Computer with INTRABEAM Software	Computer
5. Keyboard (Russian Version)	Filmy keyboard with touch panel; Keyboard (filmy)
6. Mouse	Mouse
7. Probe Adjuster Ion Chamber Holder (PAICH)	Probe centering device / holder PAICH, PAICH (probe centering device and ionization chamber holder)
8. Ionization chamber / Ionization Chamber Shield	Ionization chamber IC, Ionization chamber "IC shield"

9. Photo diode array (PDA)	Photodiode array PDA.
10. Electrometer (UNIDOS E) manufactured by PTW Freiburg Physikalisch-Technische Werkstätten Dr. Pöchlau GmbH, production site located at: Lörracher Str. 7, 79115 FREIBURG, GERMANY.	Dosimeter UNIDOS E, Electrometer UNIDOS E, UNIDOS E
11. INTRABEAM Water Phantom	Water Phantom, INTRABEAM Water Phantom
12. Isolating Transformer	---
13. INTRABEAM Workplace	INTRABEAM Workplace, INTRABEAM Cart
14. NC32 INTRABEAM Floor Stand	Floor Stand NC32 INTRABEAM, Floor Stand
15. Applicators for INTRABEAM Intraoperative Radiotherapy System in sets and on a standalone basis, with accessories P3H 2014/2108 of 01.11. 2018 r. - INTRABEAM Spherical Applicator 1.5 cm (set of INTRABEAM Spherical Applicator 1.5 cm, calibration disc) - INTRABEAM Spherical Applicator 2.0 cm (set of INTRABEAM Spherical Applicator 2.0 cm, calibration disc) - INTRABEAM Spherical Applicator 2.5 cm (set of INTRABEAM Spherical Applicator 2.5 cm, calibration disc) - INTRABEAM Spherical Applicator 3.0 cm (set of INTRABEAM Spherical Applicator 3.0 cm, calibration disc) - INTRABEAM Spherical Applicator 3.5 cm (set of INTRABEAM Spherical Applicator 3.5 cm, calibration disc) - INTRABEAM Spherical Applicator 4.0 cm (set of INTRABEAM Spherical Applicator 4.0 cm, calibration disc) - INTRABEAM Spherical Applicator 4.5 cm (set of INTRABEAM Spherical Applicator 4.5 cm, calibration disc) - INTRABEAM Spherical Applicator 5.0 cm (set of INTRABEAM Spherical Applicator 5.0 cm, calibration disc) - INTRABEAM Flat Applicator 1.0 cm (set of INTRABEAM Flat Applicator 1.0 cm, INTRABEAM Position Marker 1.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Flat Applicator 2.0 cm (set of INTRABEAM Flat Applicator 2.0 cm, INTRABEAM Position Marker 2.0cm, INTRABEAM Lumen Plug, calibration disc)	INTRABEAM Spherical Applicator, Spherical Applicator INTRABEAM Flat Applicator, Flat Applicator

- INTRABEAM Flat Applicator 3.0 cm (set of INTRABEAM Flat Applicator 3.0 cm, INTRABEAM Position Marker 3.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Flat Applicator 4.0 cm (set of INTRABEAM Flat Applicator 4.0 cm, INTRABEAM Position Marker 4.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Flat Applicator 5.0 cm (set of INTRABEAM Flat Applicator 5.0 cm, INTRABEAM Position Marker 5.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Flat Applicator 6.0 cm (set of INTRABEAM Flat Applicator 6.0 cm, INTRABEAM Position Marker 6.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Needle Applicator (set of INTRABEAM Needle Applicator, guide shaft tapered (3 pcs), guide shaft beveled (3pcs)) - INTRABEAM Surface Applicator 1.0 cm (set of INTRABEAM Surface Applicator 1.0 cm, INTRABEAM Position Marker 1.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Surface Applicator 2.0 cm (set of INTRABEAM Surface Applicator 2.0 cm, INTRABEAM Position Marker 2.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Surface Applicator 3.0 cm (set of INTRABEAM Surface Applicator 3.0 cm, INTRABEAM Position Marker 3.0cm, INTRABEAM Lumen Plug, calibration disc) - INTRABEAM Surface Applicator 4.0 cm (set of INTRABEAM Surface Applicator 4.0 cm, INTRABEAM Position Marker 4.0cm, INTRABEAM Lumen Plug, calibration disc)	INTRABEAM Needle Applicator, Needle Applicator INTRABEAM Surface Applicator, Surface Applicator
16. Container for sterilizing INTRABEAM Spherical Applicators (including re-usable filter for container)	Sterile container
17. Transport case	-
18. Dustcover	-
19. INTRABEAM DRAPE, place of production Microtek Dominicana, S.A., production site located at: Zona Franca II, I.a Romana 22000, Dominican Republic	-
20. Portable data storage device	
21. XRS Cable (cable for connection of INTRABEAM 600 and XRS radiation source)	cable XRS
22. Network/mains cable	
23. PAICH/PDA Cable	-
24. Radiation therapy treatment planning simulation system: Radiance software, manufactured by GMV Soluciones Globales	program «radiance», radiance

Internet S.A.U., production site located at: Isaac Newton 11, 28760 Tres Cantos, Madrid, Spain	
25. Operational documentation including: 25.1 INTRABEAM 600 Instruction for use 25.2 INTRABEAM Spherical Applicator Instruction for use 25.3 INTRABEAM Needle Applicator Instruction for use 25.4 INTRABEAM Surface Applicator Instruction for use 25.5 INTRABEAM Flat Applicator Instruction for use 25.6 INTRABEAM Water Phantom Instruction for use 25.7 NC32 INTRABEAM Floor Stand Instruction for use 25.8 RADIANCE INTRAOPERATIVE RADIOTHERAPY TREATMENT PLANNING Simulating SYSTEM USER's MANUAL 25.9 UNIDOS E Instruction for use	---

Annex 3. Information about indications for the use, possible side effects and contraindications.

Indications for the use:

The surgeon and/or radiation oncologist shall always carefully assess the clinical situation of the patient to determine whether the patient is suitable for IORT with INTRABEAM PRS 500 or INTRABEAM 600 system* or not

** additional information about indications for use listed in the INTRABEAM 600 Instruction for use.*

Potential side effects:

Treating patients with the radiotherapy may lead to side effects, regardless of the type of radiotherapy involved and the radiotherapy device used.

Potential side effects of the use IORT may be

- Significant wound problems such as radiation fibrosis and recurring lymphangitis
- IORT related hyperemia /erythema
- Palpable seroma
- Cellulitis

Contraindications:

There is no contraindication to IORT with INTRABEAM PRS 500 or INTRABEAM 600 system, when using INTRABEAM PRS 500 or INTRABEAM 600 system in accordance with the Operating documentation (IFU).

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Reiner Hammel


Certified Copy


I hereby certify that the above copy is a true copy of the
originals, which has been produced by me.

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