Siemens Healthcare GmbH (CT), Siemensstr. 1, 91301 Forchheim

**To the federal service for Surveillance
in Healthcare of Russian Federation**


Name	Rüdiger Körner
Department	SHS DI CT QT RCC
Telephone	+49 (9191) 18-9355
Fax	+49 (9191) 18-9782
Mobile	+49 (173) 2540949
E-mail	ruediger.koerner@siemens-healthineers.com
Date	31 August 2021

User Manual supplement for Russia

The following User Manual supplement for Russia has been printed by the manufacturer in Germany:

User Manual supplement for Russia – SOMATOM Definition Flash Computed Tomography X-ray system with accessories

Sincerely yours,
Siemens Healthcare GmbH

Rüdiger Körner
Head of Regulatory Affairs
Business Line Computed Tomography

Siemens Healthcare GmbH
Management: Bernhard Montag, President and Chief Executive Officer;
Jochen Schmitz, Christoph Zindel

Siemensstr. 1
91301 Forchheim
Germany

Tel.: +49 (9191) 18 0
siemens.com/healthcare

Chairman of the Supervisory Board: Ralf P. Thomas
Registered office: Munich, Germany; Commercial Registry: Munich, HRB 213821
WEEE-Reg.-No. DE 64872105

User manual Supplement for Russia

SOMATOM Definition Flash Computed Tomography X-ray system with accessories

1. The name of the Medical Device:

SOMATOM Definition Flash Computed Tomography X-ray system with accessories

2. Information about Manufacturer of the Medical Device

- 2.1. Designated name of the manufacturer (full name, and short name if applicable), address of the manufacturing site, post address (if differs from the official company registration):**

Legal Manufacturer: Siemens Healthcare GmbH, Henkestr. 127 DE-91052 Erlangen Germany
www.siemens.com/healthcare

Facility: Siemens Healthcare GmbH, Computed Tomography (CT), Siemensstr.1, DE-91301 Forchheim, Germany Phone: +49 9191 18-0

www.siemens.com/computedtomography

- 2.2. The name and address of the local representative in Russia:**

REPRESENTATIVE OFFICE IN RUSSIA: Siemens Healthcare LLC, Dubinskaja str. 96, 115093 Moscow, Russia, tel: +7 495 737 12 52, fax: +7 495 737 13 20.

3. Designation of the medical device: Purposes for use (Medical Device operation data)

3.1. Intended use

This computed tomography system is intended to generate and process cross-sectional images of patients by computer reconstruction of x-ray transmission data.

The images delivered by the system can be used by a trained physician as an aid in diagnosis.

3.2. Indications for Use

This computed tomography system is intended to generate and process cross-sectional images of patients by computer reconstruction of X-ray transmission data. The images delivered by the system can be used by a trained physician as an aid in diagnosis. The images delivered by the system can be used by trained staff as an aid in diagnosis, treatment preparation and radiation therapy planning. This CT system can be used for low dose lung cancer screening in high risk populations.*

*As defined by professional medical societies. Please refer to clinical literature, including the results of the National Lung Screening Trial (N Engl J Med 2011; 365:395-409) and subsequent literature, for further information.

3.3. Patient population


No patient population restrictions are known.

3.4. Contraindications

There are no known specific situations that contraindicate the use of the CT system.

3.5. Proper and improper use

Improper use means that the systems or parts of the systems are used against their intended use and indications for use.

Operation in patient vicinity

Only medical components of the system (gantry, patient table, power distribution cabinets (PDCs), accessories) are suitable for use within the patient environment.

Operation mode

The system is developed for continuous operation with intermittent loading. Thus, it should be switched off and restarted regularly. For details, see System Owner Manual.

3.6. Medical practice

The medical information presented in these instructions for use is for illustration only and is not intended to be used by the reader for the practice of medicine.

Any health care practitioner reading this medical information is reminded that they must use their own learning, training, and expertise in dealing with their individual patients. This material does not substitute for that duty and is not intended by Siemens to be used for any purpose in that regard.

The doses of the drugs mentioned are consistent with their approval labeling. The treating physician bears the sole responsibility for the diagnosis and treatment of patients, including drugs and doses prescribed in connection with such use.

3.7. Conditions for use

SOMATOM Definition Flash must only be used by persons with the necessary specialist skills who have had appropriate application training, according to country-specific regulations, for example, physicians, trained radiologists, or trained technologists.

As the user, you must have the necessary qualifications (for example, training in a relevant medical profession). You must also have been instructed in the use of the SOMATOM system (for example, detailed knowledge of the instructions for use).

SOMATOM Definition Flash can be operated in the X-ray department of a hospital or a doctor's surgery.

4. Classification of the Medical Device

4.1. Type of the MD in accordance with the nomenclature for classification:

Classification according to Annex IX, MDD

Class:	IIB
--------	-----

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According to IEC 60601

Equipment class I

Applied part
Type B (table)
Type BF (PMM)

Mode of operation
Continuous operation with intermittent loading
with duty cycle controlled by heat
unit calculators

Degree of safety
Equipment not suitable for use in the presence of a
flammable anesthetic mixture with air or with
oxygen or nitrous oxide.

According to IEC 60825

Classification of Laser Class 1M

4.2. Class of the potential risk of the medical device application:

According to GOST R 51609-2000 Medical devices.

Device is intended for transmission of energy in the form of ionizing radiation, it belong to the class 2b.

4.3. Type of the Medical Device in accordance with the Russian Nomenclature classification:

135190

4.4. Class of the Medical Device in accordance with IEC 60601-1:

In according to IEC 60601-1: 2005 class I

4.5. Type of contact with the human body:

The Medical Device is non-sterile, nonexpendable.

Device surfaces that can come into contact with patients are evaluated and proven to be biocompatible according to ISO 10993: 2009

4.6. Software classification

4.6.1 Software safety class

Class B according to IEC 62304:2006

4.6.2 Software Version

syngo CT VB10

syngo CT VB20

syngo CT VA48

4.6.3 Date of release

Version	Date
syngo CT VB10	05/2018
syngo CT VB20	05/2019
syngo CT VA48	12/2015

5. Operation (Principles of the MD operation)

A CT examination consists of the following steps:

- Patient registration
- Selecting the scan protocol
- Checking the patient position
- Gantry operator panel
- Loading the patient using the gantry control panel
- Using the Galileo 3D camera
- Setting the patient table position to zero
- Positioning the patient
- Setting the starting point for scanning
- Setting the table position to zero
- Unloading the patient
- Checking the topogram parameters
- Loading the topogram
- Acquiring the topogram
- Planning a spiral scan
- Switching on API
- Performing a spiral scan
- Reconstructing images

- Generating 3D reconstructions
- Ending the examination

Spiral scan


- Selecting the chronicle entry
- Moving the scan range
- Changing the length of the scan range
- Changing the field of view (FoV)
- Checking parameters

6. Technical Specifications of the Medical Device:

Maximum permissible deviation of technical characteristics no more than $\pm 5\%$, unless otherwise specified.

6.1. Design data:

Fig. 18: Detailed dimensions: Gantry and patient table PHS4

(A) Side view

(B) Front view

(C) Top view

(1) Slice plane

(2) Dummy plate for crane suspension

(3) 400 mm with phantom only; 350 mm with add-on components, such as head holder, table extension, etc.

(4) Movement path of tabletop support

6.2. The overview of the Medical Device:

System overview

- (1) Accessories
- (2) Scan unit
- (3) Heat exchanger
- (4) Power distribution cabinet
- (5) Image reconstruction system
- (6) Patient table
- (7) Console


Components or software functionalities that may not be part of your system configuration are nevertheless described in the present instructions for use. It may be that not all of them are marked explicitly as optional. The availability of these components or software functionalities depends on your purchase contract.


6.3. List of components and accessories of the Medical Device:


#	Name of the component in Russian	Name of the component	Description	Photo
Томограф компьютерный рентгеновский SOMATOM Definition Flash с принадлежностями: Блок базовый				
1.	Гентри (две системы трубка/детектор с приводом вращения).	Tube-Detector system consists of: -X-ray tube STRATON MX P 2 pcs -Detector 2pcs; - Slipring CPT;	- X-Ray Tube generates X-rays - Detector detects and specifies the intensity of X-ray radiation - Slipring is movable ring part of CT-system comprising scanning equipment;	

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				
6.	Блок реконструкции изображений.	IRSmx Computer for Image Reconstruction: IRSmx2d IRSmx5d CT image processor tower IRSmx5e CT image processor tower IRSmx4d IIRSmx3d	Computer for Image Reconstruction calculates and reconstructs images from the raw data obtained	 

7.	Программное обеспечение системное и прикладное: базовое программное обеспечение.	SW-Key SYNGO_General_License SW-Key Definition Flash Flash Imaging CARE Child Flash Reading FAST Spine SW-Key CARE Dose SW-Key CARE Bolus	«Flash Imaging» leading applications to boost the SOMATOM Definition Flash in redefining standards for our customers. Including SureView, Flash Spiral, Dual Adaptive Dose Shields, CARE Dose 4D, CARE kV, CARE Profile, CARE Dashboard, CARE Bolus, Dose MAP, FAST Planning, FAST Adjust, XXL Mode 2cm. CARE Child Flash Reading: Includes VRT and Workstream 4D FAST Spine SW-Key CARE Dose SW-Key CARE Bolus	
8.	Система охлаждения в составе: - внутренний блок; - наружный блок (при необходимости); - водонагреватель (при необходимости); - водопровод (при необходимости); - кабель (при необходимости).	Cooling system.	Cooling Unit is designed to maintain the necessary temperature for the smooth operation of the system.	


				
9.	Трансформатор компенсирующий.	HV transformer	High Voltage transformer feeds the system creates the necessary power	 

				
10.	Блок подключения к линии питания.	Power Distribution Cabinet: PDC A PDC B	Power Distribution Cabinet for connecting to 15 power line, supplies power to the gantry and patient table. It is used to provide power and protection.	

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



6.	Матрасы (не более 10 шт.), варианты исполнения: - стандартный; - с защитой от проливания жидкостей; - плоский; - для усиленной деки стола.	Mattress for patient table Mattress with Spill Protection patient table; Mattress with flat surface patient table Mat for High Cap.& Trauma Table Top	Mattress for patient table are used to provide easy and correct positioning of the patient on the table surface. <i>Weight: 1.8 kg</i> • Width: 470 mm • Length: 2170 mm • Height: 76 mm Mattress with improved table protection provides a secure and safe transport of the patient on the patient table. It is equipped with hinged bumpers to prevent leakage of body fluids on the table. <i>Mattress with Spill Protection Standard patient table:</i> <i>Weight: 2 kg</i> • Width: 470 mm • Length: 2170 mm • Height: 76 mm <i>Mattress with Spill Protection Multi-purpose table:</i> <i>Weight: 2,2 kg</i> <i>Dimensions:</i> • Width: 440 mm • Length: 2400 mm • Height: 76 mm Mattress with flat surface is used in clinical applications requiring a flat surface for placement of the patient. <i>Weight: 6.3 kg</i> • Width: 470 mm • Length: 2170 mm • Height: 69 mm	 
----	--	---	--	--


			Mat for high capacity & trauma table top is used for scanning nonbariatric patients on the flat, bariatric table top. Weight: 3.1 kg • Width: 495 mm • Length: 2394 mm • Height: 55 mm	 
--	--	--	---	---

7.	Подголовники/держатели для головы (не более 10 шт.).	Head holder	The head holder is used to give the correct position of the patient's head during scanning. <i>Head holder</i> Weight: 0.6 kg • Width: 226 mm • Length: 390 mm • Height: 215 mm <i>Tiltable head holder</i> Weight: 0.7 kg • Width: 220 mm • Length: 390 mm	 
8.	Опоры для головы/предплечья (не более 10 шт.).	Head-Arm-rest.	Head-Arm-Support is used to place the hands above the patient's head in a cross position. <i>Weight: 1.8 kg</i> • Width: 450 mm • Length: 480 mm • Height: 280 mm	




9.	Опоры для рук, 2 шт. (не более 10 шт.).	Arm Support, 2 pcs.	The arm supports are placed between patient and table so that the arms lie tight against the patient's body. Weight: 1.0 kg • Width: 316 mm • Length: 430 mm • Height: 150 mm	
10.	Опоры под колени (не более 10 шт.).	Knee-Leg-Support.	With the knee-leg support, you can position the lower extremities comfortably and reduce lordosis of the lumbar vertebrae. Weight: 1.2 kg • Width: 400 mm • Length: 400 mm • Height: 160 mm	

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












21.	Ключ лицензионный для активации опции коррекции артефактов от металла.	SW-Key Recon_iMAR	The iMAR CT metal artifact reduction algorithm combines three successful approaches to reduce metal artifacts: beam hardening correction (in sinogram regions of less severe metal attenuation), normalized sinogram inpainting (in sinogram regions of high metal attenuation), and frequency split (to mix back noise texture and sharp details that are potentially lost during inpainting).	
22.	Ключ лицензионный для активации опции итеративной реконструкции, варианты исполнения: - на основе сырых данных - на основе сырых данных и интеллектуальной модели	SW-Key RECON_Iterative step0 step1 Step2 Step3	Technology can reduce patient dose, and images will benefit from a higher resolution at organ borders and improved delineation of edges, e.g., to better localize lesions	

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




33.	Программно-аппаратный пакет для адаптивного объемного сканирования.	Adaptive 4D Spiral	The unique Adaptive 4D Spiral moves beyond fixed detector limitations to provide full coverage of any organ in 4D. It introduces up to 48 cm range for dynamic CTA imaging and 4D Noise Reduction to significantly improve image quality with no increase in dose or, alternately, reduce dose up to 50 % without compromising image quality.	Neurology – childhood moyamoya disease before and after surgical revascularization SOMATOM Definition Flash Collimation: 12 x 1.2 mm Adaptive 4D Spiral mode Scan time: 16 s Scan length: 140 mm Rotation time: 0.28 s 70 kV, 100 mAs CTDI_{vol}: 10.0 mAs CLP: 64.5 mAs/cm DR: 4.0 mAs Adaptive 4D Spiral delivers quantitative perfusion information of whole organs and the whole brain. Additionally the use of low kV settings will reduce the dose dramatically.  Copyright © 2014 Siemens Healthineers AG Siemens Healthineers AG Cycle time table <table><thead><tr><th>Apply</th><th>No of scans</th><th>Cycle time</th><th>Accumulated Start time</th></tr></thead><tbody><tr><td>☒</td><td>21</td><td>1.25</td><td>0.00 s</td></tr><tr><td>☒</td><td>1</td><td>1.25</td><td>26.25 s</td></tr><tr><td>☐</td><td>3</td><td>5.00</td><td></td></tr><tr><td>☐</td><td></td><td></td><td></td></tr><tr><td>☐</td><td></td><td></td><td></td></tr><tr><td>☐</td><td></td><td></td><td></td></tr><tr><td>☐</td><td></td><td></td><td></td></tr><tr><td>☐</td><td></td><td></td><td></td></tr><tr><td>☐</td><td></td><td></td><td></td></tr><tr><td>☐</td><td></td><td></td><td></td></tr><tr><td>☐</td><td></td><td></td><td></td></tr></tbody></table> No of scans: 22Exam time: 27.50 s OK Help	Apply	No of scans	Cycle time	Accumulated Start time	☒	21	1.25	0.00 s	☒	1	1.25	26.25 s	☐	3	5.00		☐				☐				☐				☐				☐				☐				☐				☐			
Apply	No of scans	Cycle time	Accumulated Start time																																																	
☒	21	1.25	0.00 s																																																	
☒	1	1.25	26.25 s																																																	
☐	3	5.00																																																		
☐																																																				
☐																																																				
☐																																																				
☐																																																				
☐																																																				
☐																																																				
☐																																																				
☐																																																				

34.	Программно-аппаратный пакет для сканирования в режиме двух энергий.	Flash Function – DE	Dual Energy scans including the Tin Filtered modes and the applications to introduce Dual Energy as part of your routine daily workflow. Includes FAST DE results and FAST DE with Workstream 4D. Item includes - syngo Dual Energy Scan with SPS - FAST DE Results #AWP - FAST DE (DE WorkStream 4D)	
35.	Программно-аппаратный пакет для сканирования сердца и оценки перфузии миокарда.	Heart Perfusion Scanning	Dynamic scan mode to visualize ventricular myocardial perfusion for identification of perfusion defects. Applies sequential technique for minimum dose and 75 ms temporal resolution for scanning even at high heart rates.	

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6.4. Technical specifications, key parameters and values, characteristics:

Maximum permissible deviation of technical characteristics no more than ± 5%, unless otherwise specified.

- Classification

Rating of the CT system according to IEC 60529 and IEC 60601

- System parameters

Rotation time at 360° turn of the gantry

- Scan mode parameters

Range for the common parameters for the typical scan modes

- Reconstruction parameters

Range for the common reconstruction parameters

- Main system components

Detailed information about the main components of the CT system

- On-site requirements

Information about the environmental requirements of the CT system

- Electromagnetic compatibility

Information about the EMC requirements of the CT system

CT system classification

Rating of the CT System according to the IEC 60529 and IEC 60601

According to:

IEC 60529

Degree of protection	IP20
----------------------	------

According to:

IEC 60601

Equipment class	I
-----------------	---

Applied part	Type B (table)
--------------	----------------

Applied part	Type BF (PMM)
--------------	---------------

Mode of operation	Continuous operation with intermittent loading
-------------------	--

Degree of safety	Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
------------------	--

6.4.1. Technical specifications, key parameters and values, characteristics VB10 version:
System parameters

Range for the common parameters for the typical scan modes.

Rotation time at 360° turn of the gantry

0.285 s	0.33 s	0.5 s	1.0 s	Tolerance ± 5%
---------	--------	-------	-------	----------------

Scan modes

Topogram

Scan range PHS 2000	20 – 2000 mm	Extended scan range
Scan time PHS 2000	0.394 – 10.76 s ²⁾	Extended scan range
Scan range MPT ²⁾	20 – 2000 mm	Extended scan range
Scan time MPT ³⁾	0.394 – 10.76 s ⁴⁾	Extended scan range
Views	a.p., p.a., lateral	

Sequence

General

Number of uninterrupted scans per range	100	
Tomogram scan time		Tolerance ± 10%
Temporal resolution with dedicated reconstruction	75 ms	0.28 s/Dual tube
i-Sequence (optional)		
Scan time		max. 100 s

Spiral

General

Scan time	max. 160 s ^{a)}	Tolerance: ± 10%
Scan length	max. 1840 mm	
Delay time	2 - 600 s	
Temporal resolution	75 ms	biphasic/0.28
Max. number of auto-ranges per protocol	33	

^{a)} optional

- 1) Tolerance ± 10%
- 2) Optional
- 3) Optional
- 4) Tolerance ± 10%

Spiral	Pitch	Pitch increment
Routine	0.35 – 1.5	0.05
Flash (Dual source high pitch spiral)	1.55 – 3.2	0.05
Neuro	0.35 – 1.5	0.05
Obese	0.35 – 1.5	0.05
UHR / z-UHR ^{a)}	0.3 – 1.0	0.05
Chest pain	0.17 – 1.5	0.01
Chest Pain Flash	1.55 – 3.2	0.05
Cardio ^{a)}	0.17 – 1.5	0.01
Cardio Flash ^{a)}	2.8 – 3.4	0.05
Respiratory ^{a)}	0.09 – 0.2	0.01
i-Spiral ^{a)}	0.5 – 1.5	0.05
Dual Energy ^{a)}	0.2 – 1.2	0.05
Dual Energy Cardio ^{a)}	0.17 – 1.5	0.01

^{a)} Optional

Multiscan

Scan time	max. 100 s	Tolerance: ± 10%
-----------	------------	------------------

Adaptive 4D spiral (optional)

Adaptive 4D spiral	For cycle time	Scan length
Perfusion	1.5 s	100 mm
Multiphase	1.5 s	150 mm
	2.0 s	216 mm
	2.5 s	284 mm
Fast CTA	2.5 s	381 mm
	3.0 s	481 mm
Topogram scan time $\pm$ 10%		

Scan length: Nominal range to be displayed in the range UI item.
For cycle time: Nominal time to be displayed in the range UI item.

Imaging for RTP

Test conditions

Unless specified otherwise in the individual methods, the tests have been conducted under the following ambient / operating conditions specified by the manufacturer.

Ambient temperature	Relative humidity	Barometric pressure
18 - 28 °C	20 - 75 %	800 - 1060 hPa

Patient table sag test according to IEC 60601 2-44

According to IEC 60601 2 44, Clause 201.101.4

Test is conducted according to IEC 60601-2-44 Clause 201.101.4. For that purpose, a load of 135 kg was distributed over a length of 1.9 m. The table height in the scan plane is measured at the gantry-side end of the table top (i.e. 0 mm) and in steps of 400 mm. The deviation is stated as height difference between the adjacent measuring positions.

The actual values may vary depending on load, load distribution and table top or overlay. It is recommended to perform analogous measurements with the case specific setup if table sag corrections are considered in the clinical workflow of Radiation Therapy Planning.

Multi Purpose Table with RTP table top or compatible RTP overlays

Measurement position	Height deviation to previous position
0 mm	initial reference position
400 mm	< 1 mm
800 mm	< 1 mm
1200 mm	< 2 mm

Patient table 1600 with 3rd party RTP overlays*

Measurement position	Height deviation to previous position
0 mm	initial reference position
400 mm	< 5 mm
800 mm	< 2 mm
1200 mm	< 3 mm

* ask your sales contact for a current list of compatible overlays

Multi-purpose table Sliding Gantry with RTP table top

Measurement position	Height deviation to previous position
0 mm	initial reference position
400 mm	< 2 mm
800 mm	< 3 mm
1200 mm	< 3 mm

Patient table sag test (AAPM Radiation Therapy Committee)

According to the Report of the AAPM Radiation Therapy Committee Task Group
No. 66, APPENDIX E

The table top is loaded with 75kg distributed to simulate the weight of a patient.
 The table top height in the scan plane is measured at the gantry-side end of the table top. The height difference of the table top surface in the scan plane is below 2 mm for horizontal table movement within the movement ranges, according to the tables below.
 The actual values may vary depending on load, load distribution and table top or overlay.

Table type	Movement range	Height deviation
Multi-purpose table	1400 mm	≤ 2 mm

The actual values may vary depending on load, load distribution and table top or overlay.

Table type	Movement range	Height deviation
Patient table 1600	560 mm	≤ 2 mm

Light markers

According to IEC 60601-2-44, Clause 201.101.5
 The integral light markers fulfill IEC 60601-2-44, Clause 201.101.5 and can be used for patient marking in RTP.

Geometric accuracy of image data

Angular Alignment of CT Images
 Test method according to IEC 60601-2-44, Subclause 201.101.8.3.
 The deviation over a distance of 200 mm is less than 3 mm.

Accuracy of z-position for helical scans
 Test method according to IEC 60601-2-44, Subclause 201.101.8.4.
 The deviation over a distance of 300 mm is less than 1mm.

System reconstruction parameters

Range for the common parameters for the typical scan modes.

Image reconstruction

Slice width	0.5 – 10 mm	
Slice increment	0.1 – 10 mm	
Scan field	50 cm	For high pitch modes 33 cm
Recon field (body)	5 – 50 cm	
Recon field (head)	5 – 30 cm	
Recon matrix	512 x 512 pixel	
HU scale	-1024 to + 3071	1 HU increment
Extended HU scale	-10240 to +30710	10 HU increments

Image reconstruction (optional)

Extended recon field (adult body)	50 – 78 cm	Extended field of view
Slice width	0.5 – 6 mm	UHR
Slice width	0.4 – 6 mm	2-UHR

Main system components


Detailed information about the main components of the CT system.

Tube and generator

2x X-ray tube

For X-ray properties

Cooling curve for the tube housing assembly

x Time [s]
y Rel. heat content [%]


The tube housing assembly with its closed-loop cooling system integrated in the gantry has a continuous heat transfer capability of up to 4.6 kW.

Due to the anode cooling concept, the heat storage capacity is not limited by the anode heat storage capacity, but by the cooling system.

Tube rating charts

The following graph shows the tube rating chart for focal spot UHR 0,7 x 0,7 mm* and tube rating chart for flying focal spots 0,9 x 1,1 mm*.

*) According to: IEC 60336

- x Time [s]
y Scan power [kW]
1 Cooling system limit
2 load limit 120 Hz
3 load limit 140 Hz
4 load limit rpp 1x1, 160 Hz
5 load limit rpp 2x1, 160 Hz
6 load limit rpp 2x2, 160 Hz
7 load limit rpp 1x2, 160 Hz

Output current

Voltage [kV]	current range [mA]
70	60 ... 500
80	20 ... 650
100	20 ... 650
120	20 ... 800
140	20 ... 714

Tolerance tube voltage $\pm 10\%$

Tolerance tube current $\pm 20\%$

Detector

Detector (360° scan)

Number of projections	1152, 2304, 4608
Elements (tube A)	47104 (736 x 64)
Elements (tube B)	30720 (480 x 64)
Number of detector rows	2 x 64

Gantry

Aperture	78 cm	
Distance gantry front to scan plane	35 cm	
Distance focal spot to isocenter	595 mm	
Scan field diameter (tube A)	50 cm	
Scan field diameter (tube B)	35 cm	
Scan plane lightmarker	accuracy lightmarker to scan plane	Tolerance: ± 1 mm
Gantry front lightmarker *	accuracy lightmarker to scan plane	Tolerance: ± 1 mm
Sagittal plane lightmarker	accuracy lightmarker to vertical plane through isocenter	Tolerance: ± 1 mm

***) optional**

Coronal plane lightmarker	accuracy lightmarker to horizontal plane through isocenter	Tolerance: ± 2 mm
---------------------------	--	-------------------

***) optional**

Patient table 2000

Vertical patient load position ⁵⁾	490 mm	Lower bound tolerance: - 5 mm Upper bound tolerance: + 5 mm
Max. table load	227 kg / 500 lbs	
Max. table movement after emergency stop (during patient positioning)		≤ 10 mm

Multi purpose table (optional)

Vertical patient load position ⁶⁾	550 mm	Lower bound tolerance: - 5 mm Upper bound tolerance: + 5 mm
Max. table load	227 kg / 500 lbs	
Max. table load with bariatric/trauma table top	307 kg / 676 lbs	
Max. table movement after emergency stop (during patient positioning)		≤ 10 mm

Image display

Monitor resolution	1280 × 1024 pixel
Image display matrix	max. 1024 × 1024 pixel

Storage

	ICS main storage	approx. 520000 images
IRS main storage		
IRSmx4d	x	
3800 GB		
IRSmx5d	x ^{o)}	
4000 GB		
IRSmx5e	x ^{o)}	
4000 GB		
^{o)} optional		

- 5) At table top
6) At table top

CD Media	700 MB (1200 images)
DVD Media	4.7 GB (8000 images)

Laser localizer

Classification in accordance with IEC 60825	Class 1M
Spectral range	635nm ±10nm
Output power	≤ 3mW

On-site requirements

System connected values

Line voltage and line current

Line voltage	Maximum line current (peak)
380 V ± 10%	230 A per phase
400 V ± 10%	225 A per phase
420 V ± 10%	220 A per phase
440 V ± 10%	215 A per phase
460 V ± 10%	210 A per phase
480 V ± 10%	205 A per phase

Line voltage, nominal line frequency

Nominal Voltage 3/N-	380-480 V in 20 V steps
Nominal line frequency	50; 60 Hz

Technique factors for max. line current

Tube voltage	140 kV
Tube current	714 mA
Scan time	max. 4 s

Power consumption

Operating mode	Power consumption	
	Standard CT system	With optional water/air split cooling system
Operation / Input power	283 kVA (140 kV, 714 mA, 4 s, dual tube)	293 kVA (140 kV, 714 mA, 4 s, dual tube)
System ON (stand-by)	6 kVA	16 kVA
System ON (Computer On)	3 kVA	3 kVA
System OFF	1 kVA	1 kVA

Environment requirements

Climatic conditions inside and outside of the examination room

Climatic conditions inside and outside of the examination room

Relative humidity	20 – 75 %	
Temperature gradient	max. 6 K/h	operation
Atmospheric pressure	800 hPa – 1060 hPa	
Temperature range	18 – 28 °C (64.4 – 82.4 °F)	

CT system and its components inside and outside of the examination room

6.4.2. Technical specifications, key parameters and values, characteristics VB20 version: System parameters

Range for the common parameters for the typical scan modes.

Rotation time at 360° turn of the gantry

0.285 s	0.33 s	0.5 s	1.0 s	Tolerance ± 5%
---------	--------	-------	-------	----------------

Scan modes

Topogram

Scan range PHS 2000	20 – 2000 mm	Extended scan range
Scan time PHS 2000	0.394 – 10.76 s ¹⁾	Extended scan range
Scan range MPT ²⁾	20 – 2000 mm	Extended scan range
Scan time MPT ²⁾	0.394 – 10.76 s ¹⁾	Extended scan range
Views	a.p., p.o., lateral	

Sequence

General

Number of uninterrupted scans per range	100	
Tomogram scan time	Tolerance: ± 10%	
Temporal resolution with dedicated reconstruction	75 ms	0.28 s ¹⁾

¹⁾Dual tube

i-Sequence (optional)

Scan time	max. 100 s
-----------	------------

Spiral

General

Scan time	max. 160 s ¹⁾	Tolerance: ± 10%
Scan length	max. 1840 mm	
Delay time	2 - 600 s	
Temporal resolution	75 ms	biphasic/0.28
Max. number of auto-ranges per protocol	17	

1) Tolerance ± 10%

2) Optional

Spiral	Pitch	Pitch increment
Routine	0.35 – 1.5	0.05
Flash (Dual source high pitch spiral)	1.55 – 3.2	0.05
Neuro	0.35 – 1.5	0.05
Obese	0.35 – 1.5	0.05
UHR / z-UHR ¹⁾	0.3 – 1.0	0.05
Chest pain	0.17 – 1.5	0.01
Chest Pain Flash	1.55 – 3.2	0.05
Cardio ⁴	0.17 – 1.5	0.01
Cardio Flash ¹⁾	2.8 – 3.4	0.05
Respiratory ²⁾	0.09 – 0.2	0.01
i-Spiral ³⁾	0.5 – 1.5	0.05
Dual Energy ⁴	0.2 – 1.2	0.05
Dual Energy Cardio ¹⁾	0.17 – 1.5	0.01

¹⁾ Optional

Multiscan

Scan time	max. 100 s	Tolerance: ± 10%
-----------	------------	------------------

Adaptive 4D spiral (optional)

	For cycle time ¹⁾	Scan length ²⁾
Perfusion	1.5 s	100 mm
	1.5 s	150 mm
Multiphase	2.0 s	216 mm
	2.5 s	284 mm
Fast CTA	2.5 s	381 mm
	3.0 s	481 mm

Topogram scan time ± 10%

¹⁾ Nominal time to be displayed in the range UI item; ²⁾ Nominal range to be displayed in the range UI item

Scan length: Nominal range to be displayed in the range UI item.

For cycle time: Nominal time to be displayed in the range UI item.

Imaging for RTP

Test conditions

Unless specified otherwise in the individual methods, the tests have been conducted under the following ambient / operating conditions specified by the manufacturer.

Ambient temperature	Relative humidity	Barometric pressure
18 - 28 °C	20 - 75 %	800 - 1060 hPa

Patient table sag test (IEC 60601-2-44)

According to IEC 60601 2 44, Clause 201.101.4

Test is conducted according to IEC 60601-2-44 Clause 201.101.4. For that purpose, a load of 135 kg was distributed over a length of 1.9 m. The table height in the scan plane is measured at the gantry-side end of the table top (i.e. 0 mm) and in steps of 400 mm. The deviation is stated as height difference between the adjacent measuring positions.

The actual values may vary depending on load, load distribution and table top or overlay. It is recommended to perform analogous measurements with the case specific setup if table sag corrections are considered in the clinical workflow of Radiation Therapy Planning.

Multi Purpose Table with RTP table top or compatible RTP overlays

Measurement position	Height deviation to previous position
0 mm	initial reference position
400 mm	< 1 mm
800 mm	< 1 mm
1200 mm	< 2 mm

Patient table sag test (AAPM Radiation Therapy Committee)

According to the Report of the AAPM Radiation Therapy Committee Task Group No. 66, APPENDIX E

The table top is loaded with 75kg distributed to simulate the weight of a patient.

The table top height in the scan plane is measured at the gantry-side end of the table top. The height difference of the table top surface in the scan plane is below 2 mm for horizontal table movement within the movement ranges, according to the tables below.

The actual values may vary depending on load, load distribution and table top or overlay.

Table type	Movement range	Height deviation
Multi-purpose table	1400 mm	≤ 2 mm

Light markers

According to IEC 60601-2-44, Clause 201.101.5

The integral light markers fulfill IEC 60601-2-44, Clause 201.101.5 and can be used for patient marking in RTP.

Geometric accuracy of image data

Angular Alignment of CT Images

Test method according to IEC 60601-2-44, Subclause 201.101.8.3.

The deviation over a distance of 200 mm is less than 3 mm.

Accuracy of z-position for helical scans

Test method according to IEC 60601-2-44, Subclause 201.101.8.4.

The deviation over a distance of 300 mm is less than 1mm.

System reconstruction parameters

Range for the common parameters for the typical scan modes.

Image reconstruction

Slice width	0.5 – 10 mm	
Slice increment	0.1 – 10 mm	
Scan field	50 cm	For high pitch modes 33 cm
Recon field (body)	5 – 50 cm	
Recon field (head)	5 – 30 cm	
Recon matrix	512 x 512 pixel	
HU scale	-1024 to +3071	1 HU increment
Extended HU scale	-10240 to +30710	10 HU increments

Image reconstruction (optional)

Extended recon field (adult body)	50 – 78 cm	Extended field of view
Slice width	0.5 – 6 mm	UHR
Slice width	0.4 – 6 mm	z-UHR

Main system components


Detailed information about the main components of the CT system.

Tube and generator

2x X-ray tube

For X-ray properties

Cooling curve for the tube housing assembly

x Time [s]
y Rel. heat content (%)


The tube housing assembly with its closed-loop cooling system integrated in the gantry has a continuous heat transfer capability of up to 4.6 kW.

Due to the anode cooling concept, the heat storage capacity is not limited by the anode heat storage capacity, but by the cooling system.

Tube rating charts

The following graph shows the tube rating chart for focal spot UHR 0,7 x 0,7 mm* and tube rating chart for flying focal spots 0,9 x 1,1 mm*.

*) According to: IEC 60336

x Time [s]
y Scan power [kW]
1 Cooling system limit

- 2 load limit 120 Hz
- 3 load limit 140 Hz
- 4 load limit rpp 1x1, 160 Hz
- 5 load limit rpp 2x1, 160 Hz
- 6 load limit rpp 2x2, 160 Hz
- 7 load limit rpp 1x2, 160 Hz

Output current

Voltage [kV] ¹⁾	Current range [mA] ²⁾
70	60 ... 500
80	20 ... 650
100	20 ... 650
120	20 ... 800
140	20 ... 714

¹⁾ Tolerance tube voltage $\pm 10\%$; ²⁾ Tolerance tube current $\pm 20\%$

Detector

Detector (360° scan)	
Number of projections	1152, 2304, 4608
Elements (tube A)	47104 (736 x 64)
Elements (tube B)	30720 (480 x 64)
Number of detector rows	2 x 64

Gantry

Aperture	78 cm
Distance gantry front to scan plane	35 cm
Distance focal spot to isocenter	595 mm
Scan field diameter (tube A)	50 cm
Scan field diameter (tube B)	35 cm
Scan plane lightmarker	accuracy lightmarker to scan plane ²⁾
Gantry front lightmarker ¹⁾	accuracy lightmarker to scan plane ²⁾
Sagittal plane lightmarker	accuracy lightmarker to vertical plane through isocenter ²⁾
Coronal plane lightmarker	accuracy lightmarker to horizontal plane through isocenter ³⁾

¹⁾ Optional; ²⁾ Tolerance: $\pm 1.0^\circ$; ³⁾ Tolerance: $\pm 0.5^\circ$

Patient table 2000

Vertical patient load position ¹⁾	480 mm ²⁾
Max. table load	227 kg / 500 lbs
Max. table movement after emergency stop (during patient positioning)	≤ 10 mm
Max. table movement after emergency stop (during scan modes)	≤ 25 mm
Max. table movement after power failure	≤ 50 mm

¹⁾ At table top; ²⁾ Lower bound tolerance: -5 mm; Upper bound tolerance: +5 mm

Multi purpose table (optional)

Vertical patient load position ¹⁾	550 mm ²⁾
Max. table load	227 kg / 500 lbs
Max. table load with bariatric/trauma table top	307 kg / 676 lbs
Max. table movement after emergency stop (during patient positioning)	≤ 10 mm
Max. table movement after emergency stop (during scan modes)	≤ 25 mm
Max. table movement after power failure	≤ 50 mm

¹⁾ At table top; ²⁾ Lower bound tolerance: -5 mm; Upper bound tolerance: +5 mm

Image display

Monitor resolution	1280 × 1024 pixel
Image display matrix	512 × 512 pixel

Storage

	ICS Tower 12	ICS Tower 13
ICS main storage	approx. 550,000 images	approx. 2.1 million images

¹⁾ At standard matrix size of 512 × 512 pixel

	Main storage
IRSmx4d	3,800 GB
IRSmx5e ¹⁾	4,000 GB
CD Media	700 MB (1,200 images)
DVD Media	4.7 GB (8,000 images)

¹⁾ Optional

Laser localizer

Classification in accordance with IEC 60825	Class 1M
Spectral range	635nm ±10nm

Output power	≤ 3mW
--------------	-------

On-site requirements

System connected values

Line voltage and line current

Line voltage	Maximum line current (peak)
380 V ± 10%	230 A per phase
400 V ± 10%	225 A per phase
420 V ± 10%	220 A per phase
440 V ± 10%	215 A per phase
460 V ± 10%	210 A per phase
480 V ± 10%	205 A per phase

Nominal voltage and nominal line frequency

Nominal Voltage 3/N~	380–480 V in 20 V steps
Nominal line frequency	50; 60 Hz

Technique factors for max. line current

Tube voltage	140 kV
Tube current	714 mA
Scan time	max. 4 s

Power consumption

Operating mode	Power consumption	
	Standard CT system	With optional water/air split cooling system
Operation / Input power dual tube	283 kVA (140 kV, 714 mA, 4 s)	293 kVA (140 kV, 714 mA, 4 s)
System ON (stand-by)	6 kVA	16 kVA
System ON (Computer On)	3 kVA	3 kVA
System OFF	1 kVA	1 kVA

Highest power output of high-voltage generator (for 4 s)

Operating mode	Power output
Operation / Input power single tube	96 kW (120 kV, 800 mA, 4 s)
Operation / Input power dual tube	192 kW (120 kV, 800 mA, 4 s)

European Coordination Committee of the Radiological, Electromedical and Healthcare IT Industry (COCIR) describes a 24-h power consumption model for CT systems. The Green Public Procurement (GPP), on the other hand, describes a simplified version of the COCIR consumption model.

Power consumption within 24 h

24-h use scenario	24-h power consumption ¹⁾	
	COCIR model	GPP model
Off scenario ²⁾	42.9 kWh	43.1 kWh
Low power scenario	49.8 kWh	49.9 kWh
Idle scenario	83.1 kWh	83.2 kWh

¹⁾ Values may vary approx. +/- 3% due to specific system conditions, for example charging of UPS, etc.; ²⁾ With Site On/Off switch

Environment requirements

Climatic conditions inside and outside of the examination room

Relative humidity	20 – 75 %
Temperature gradient	max. 6 K/h (operation)
Atmospheric pressure	800 – 1060 hPa
Temperature range	18 – 28 °C (64.4 – 82.4 °F)

6.4.3. Technical specifications, key parameters and values, characteristics VA48 version: System parameters

Range for the common parameters for the typical scan modes.

Rotation time at 360° turn of the gantry

0.28 s	0.33 s	0.5 s	1.0 s	Tolerance $\pm 5\%$
--------	--------	-------	-------	---------------------

Scan modes

Topogram

Scan range PHS 2000	128 – 1970 mm	Extended scan range
Scan time PHS 2000	1.5 – 20.0 s ¹⁾	Extended scan range
Scan range MPT ²⁾	128 – 1970 mm	Extended scan range
Scan time MPT ³⁾	1.5 – 20.0 s ⁴⁾	Extended scan range
Views	a.p., p.a., lateral	

Sequence

General

Number of uninterrupted scans per range	100	
Tomogram scan time		Tolerance $\pm 10\%$
Temporal resolution with dedicated reconstruction	75 ms	0.28 s/Dual tube - Quickscan

i-Sequence (optional)

Scan time	max. 100 s
-----------	------------

Spiral

General

Scan time	max. 160 s*	Tolerance: $\pm 10\%$
Scan length	max. 1840 mm	
Delay time	2 - 600 s	
Temporal resolution	75 ms	biphasic/0.28
Max. number of auto-ranges per protocol	33	

*) optional

- 1) Tolerance $\pm 10\%$
- 2) Optional
- 3) Optional
- 4) Tolerance $\pm 10\%$

Spiral	Pitch	Pitch increment
Routine	0.35 – 1.5	0.05
Flash (Dual source high pitch spiral)	1.55 – 3.2	0.05
Neuro	0.35 – 1.5	0.05
Obese	0.35 – 1.5	0.05
UHR / z-UHR*	0.3 – 1.0	0.05
Chest pain	0.17 – 1.5	0.01
Chest Pain Flash	1.55 – 3.2	0.05
Cardio*	0.17 – 1.5	0.01
Cardio Flash*	2.8 – 3.4	0.05
Respiratory*	0.09 – 0.2	0.01
i-Spiral*	0.5 – 1.5	0.05
Dual Energy*	0.2 – 1.2	0.05
Dual Energy Cardio*	0.17 – 1.5	0.01
*) Optional		

Multiscan

Scan time	max. 100 s	Tolerance: $\pm 10\%$
-----------	------------	-----------------------

Adaptive 4D spiral (optional)

Adaptive 4D spiral	For cycle time	Scan length
Perfusion	1.0 s	70 mm
	1.25 s	84 mm
	1.5 s	100 mm
Multiphase	1.5 s	150 mm
	1.75 s	182 mm
	2.0 s	216 mm
	2.25 s	251 mm
	2.5 s	284 mm
Fast CTA	2.5 s	381 mm
	3.0 s	481 mm
Topogram scan time $\pm 10\%$		

Scan length: Nominal range to be displayed in the range UI item.

For cycle time: Nominal time to be displayed in the range UI item.

Imaging for RTP

Test conditions and references

Unless specified otherwise in the individual methods, the tests have been conducted under the following ambient / operating conditions specified by the manufacturer.

Ambient temperature	Relative humidity	Barometric pressure
18 - 28 °C	20 - 75 %	800 - 1060 hPa

Patient table sag test (IEC 60601-2-44)

According to IEC 60601 2 44, Clause 201.101.4

Test is conducted according to IEC 60601-2-44 Clause 201.101.4. For that purpose, a load of 135 kg was distributed over a length of 1.9 m. The table height in the scan plane is measured at the gantry-side end of the table top (i.e. 0 mm) and in steps of 400 mm. The deviation is stated as height difference between the adjacent measuring positions.

The actual values may vary depending on load, load distribution and table top or overlay. It is recommended to perform analogous measurements with the case specific setup if table sag corrections are considered in the clinical workflow of Radiation Therapy Planning.

Multi Purpose Table with RTP table top or compatible RTP overlays

Measurement Position	Height Deviation to previous position
0 mm	n.a. - initial reference position
400 mm	< 1 mm
800 mm	< 1 mm
1200 mm	< 2 mm

Patient table sag test (AAPM Radiation Therapy Committee)

According to the Report of the AAPM Radiation Therapy Committee Task Group No. 66, APPENDIX E

The table top is loaded with 75kg distributed to simulate the weight of a patient.

The table height in the scan plane is measured at the gantry-side end of the table top. The height difference of the table top surface in the scan plane is below 2mm for horizontal table movement.

Table Type	Travel Range	Height Deviation
Multi Purpose Table	1400 mm	≤ 2 mm

The actual values may vary depending on load, load distribution and table top or overlay.

Light markers

According to IEC 60601-2-44, Clause 201.101.5

The integral light markers fulfill IEC 60601-2-44, Clause 201.101.5 and can be used for patient marking in RTP.

Geometric accuracy of image data

Angular Alignment of CT Images

Test method according to IEC 60601-2-44, Subclause 201.101.8.3.

The deviation over a distance of 200 mm is less than 3 mm.

Accuracy of z-position for helical scans

Test method according to IEC 60601-2-44, Subclause 201.101.8.4.

The deviation over a distance of 300 mm is less than 1mm.

System reconstruction parameters

Range for the common parameters for the typical scan modes.

Image reconstruction

Slice width	0.5 – 10 mm	
Slice increment	0.1 – 10 mm	
Scan field	50 cm	For high pitch modes 33 cm
Recon field (body)	5 – 50 cm	
Recon field (head)	5 – 30 cm	
Recon matrix	512 x 512 pixel	
HU scale	-1024 to + 3071	1 HU increment
Extended HU scale	-10240 to +30710	10 HU increments

Image reconstruction (optional)

Extended recon field (adult body)	50 – 78 cm	Extended field of view
Slice width	0.5 – 6 mm	UHR
Slice width	0.4 – 6 mm	z-UHR

Main system components


Detailed information about the main components of the CT system.

Tube and generator

2x Straton MX P

For X-ray properties

Cooling curve for the tube housing assembly

x Time [s]
y Rel. heat content [%]


The tube housing assembly with its closed-loop cooling system integrated in the gantry has a continuous heat transfer capability of up to 4.6 kW.

Due to the anode cooling concept, the heat storage capacity is not limited by the anode heat storage capacity, but by the cooling system.

Tube rating charts

The following graph shows the tube rating chart for focal spot UHR 0,7 x 0,7 mm* and tube rating chart for flying focal spots 0,9 x 1,1 mm*.

*) According to: IEC 60336

- x Time [s]
y Scan power [kW]
1 Cooling system limit
2 load limit 120 Hz
3 load limit 140 Hz
4 load limit rpp 1x1, 160 Hz
5 load limit rpp 2x1, 160 Hz
6 load limit rpp 2x2, 160 Hz
7 load limit rpp 1x2, 160 Hz

Output current

Voltage [kV]	70	80	100	120	140
min. [mA]	60	20	20	20	20
max. 100 kW [mA]	500	650	650	800	714

Tolerance tube voltage $\pm 10\%$

Tolerance tube current $\pm 20\%$

Detector

Detector (360° scan)

Number of projections	1152, 2304, 4608
Elements (tube A)	47104 (736 x 64)
Elements (tube B)	30720 (480 x 64)
Number of detector rows	2 x 64

Gantry

Aperture	78 cm	
Distance gantry front to scan plane	35 cm	
Distance focal spot to isocenter	595 mm	
Scan field diameter (tube A)	50 cm	
Scan field diameter (tube B)	35 cm	
Scan plane lightmarker	accuracy lightmarker to scan plane	Tolerance: ± 1 mm
Gantry front lightmarker	accuracy lightmarker to scan plane	Tolerance: ± 1 mm
Sagittal plane lightmarker	accuracy lightmarker to vertical plane through isocenter	Tolerance: ± 2 mm
Coronal plane lightmarker	accuracy lightmarker to horizontal plane through isocenter	Tolerance: ± 2 mm

Patient table 2000

Vertical patient load position ⁵⁾	490 mm	Lower bound tolerance: - 5 mm Upper bound tolerance: + 5 mm
Max. table load	227 kg / 500 lbs	
Max. table movement after emergency stop (during patient positioning)		≤ 10 mm

Multi purpose table (optional)

Vertical patient load position ⁶⁾	550 mm	Lower bound tolerance: - 5 mm Upper bound tolerance: + 5 mm
Max. table load	227 kg / 500 lbs	
Max. table load with bariatric/trauma table top	307 kg / 676 lbs	
Max. table movement after emergency stop (during patient positioning)		≤ 10 mm

5) At table top

6) At table top

Image display

Monitor resolution	1280 × 1024 pixel
Image display matrix	max. 1024 × 1024 pixel

Storage

ICS main storage	approx. 520000 images
IRS main storage	min. 3.8 TB
CD Media	700 MB (1200 images)
DVD Media	4.7 GB (8000 images)

Laser localizer

Classification in accordance with IEC 60825	Class 1M
Spectral range	635nm $\pm$ 10nm
Output power	$\leq$ 3mW

On-site requirements

System connected values

Line voltage and line current

Line voltage	Maximum line current (peak)
380 V $\pm$ 10%	230 A per phase
400 V $\pm$ 10%	225 A per phase
420 V $\pm$ 10%	220 A per phase
440 V $\pm$ 10%	215 A per phase
460 V $\pm$ 10%	210 A per phase
480 V $\pm$ 10%	205 A per phase

Line voltage, nominal line frequency

Nominal Voltage 3/N~	380–480 V in 20 V steps
Nominal line frequency	50; 60 Hz

Technique factors for max. line current

Tube voltage	140 kV
Tube current	714 mA
Scan time	max. 4 s

Power consumption

Operating mode	Power consumption	
	Standard CT system	With optional water/air split cooling system
Operation / Input power	283 kVA (140 kV, 714 mA, 4 s, dual tube)	293 kVA (140 kV, 714 mA, 4 s, dual tube)
System ON (stand-by)	6 kVA	16 kVA
System ON (Computer On)	3 kVA	3 kVA
System OFF	1 kVA	1 kVA

Environment requirements

Climatic conditions inside and outside of the examination room

Relative humidity	20 – 75 %	
Temperature gradient	max. 6 K/h	operation
Atmospheric pressure	800 hPa – 1060 hPa	
Temperature range	18 – 28 °C (64.4 – 82.4 °F)	


CT system and its components inside and outside of the examination room

6.5 Labeling Information and requirements

6.5.1 Labeling Information and requirements VB10 version

The numbers and parameters on the labels are examples only. Always refer to the original labels that are applied on the CT system and system components.

(1) Approval label (Saudi Arabia only, MOH)

(2) Follow the Instructions for Use

(1)

Approval label (Saudi Arabia only, MOH)

MOH – Ministry of Health


(2)

Follow the Instructions for Use

The user must consult and follow the Instructions for Use before using the CT system and its equipment.
An electronic copy of the Instructions for Use is available at <http://www.siemens.com/medImaging-Manuals>.

Rear side

- (1) Warning, check for absence of power
- (2) Product identification label (System)
- (3) Product identification label (Gantry)
- (4) Approval label (USA only, DHHS)
- (5) Approval label (USA and Canada, NRTL)
- (6) Radio frequency information (Only for systems with an optional wireless i-Control)

- (1)
Warning, check for absence of power

Warning

Even when the CT system is switched off there is voltage present. To ensure safe working conditions, all external live lines must be switched off on-site and secured against being turned on again. Before beginning work wait 2 minutes and check for absence of power.

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(4)

Approval label (USA only, DHHS)

This product complies with DHHS regulations
21 CFR Subchapter J, applicable at date of
manufacture.

Manufactured: XXXXXXXXXXXX

Siemens Healthcare GmbH
Henkestr. 127, 91052 Erlangen
Germany

DHHS – Department of Health and Human Services

(5)


Approval label (USA and Canada, NRTL)

NRTL – Nationally Recognized Testing Laboratory

(6)

Radio frequency information


Only for systems with an optional wireless i-Control.

The individual, country specific approval may vary on this label.

Gantry frame

The labels at the gantry frame are located beneath the gantry cover.

- (1) Product status label
- (2) Product identification label (Gantry)
- (3) Approval label (USA only, DHHS)**
- (4) Approval label (USA and Canada, NRTL)**
- (5) Product identification label (System)**

(1)


Product status label

Rev.	00	01	02	03	04	05	06	07	08	09	10	11
12	13	14	15	16	17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32	33	34	35	36	37

Refer to the original label that is applied on the device for information about the current product version. This label is an example only.


(2)

Product identification label

(3)

Approval label (USA only, DHHS)

DHHS – Department of Health and Human Services

(4)


Approval label (USA and Canada, NRTL)


NRTL – Nationally Recognized Testing Laboratory

(5)

Product identification label (System)

Tube housing

- (1) Product status label
- (2) Focal spot size
- (3) Product identification label (Tube assembly)
- (4) Approval label (USA, DHHS)
- (5) Approval label (USA and Canada, NRTL)
- (6) Product identification label (X-ray tube)

(1)


Product status label

Rev.	00	01	02	03	04	05	06	07	08	09	10	11
12	13	14	15	16	17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32	33	34	35	36	37

Refer to the original label that is applied on the device for information about the current product version. This label is an example only.

(2)

Focal spot size

(3)

Product identification label


– or –

Product identification label

(4)


Approval label (USA, DHHS)

DHHS – Department of Health and Human Services

– or –


Approval label (USA, DHHS)

DHHS – Department of Health and Human Services

(5)

Approval label (USA and Canada, NRTL)

NRTL – Nationally Recognized Testing Laboratory

(6)

Product identification label

SIEMENS	
(1P) Model No.	XXXXXXX
(S) Serial No.	XXXXXXXXXX
(1P) Model No.	
(S) Serial No.	
Roehre/Tube	STRATON XX X
Siemens Aktiengesellschaft, Wittelsbacher Platz2, D-8033 Muenchen (Germany)	
Made in Germany	


– or –

Product identification label

SIEMENS	
STRATON XX	
	(01) xxxxxxxxxxxxxxxx GTIN
	(240) xxxxxxxx MODEL
	(422) xxxxx ORIGIN DE
	(21) xxxxx SERIAL
Siemens Healthcare GmbH Henkestr. 127, 91052 Erlangen, Germany	
Made in Germany	

Tube collimator

Handwritten signature

- (1) RoHS label
- (2) Approval label (USA, DHHS)
- (3) Filtering label
- (4) Product identification label
- (5) Product status label

(1)
RoHS label

(2)
Approval label (USA, DHHS)

**This product complies with DHHS regulations
21 CFR Subchapter J, applicable at date of
manufacture.**

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DHHS – Department of Health and Human Services

– or –

Approval label (USA, DHHS)

This product complies with DHHS regulations
21 CFR Subchapter J, applicable at date of
manufacture.

Manufactured: XXXXXXXXXXXXX

Siemens Healthcare GmbH
Henkestr. 127, 91052 Erlangen
Germany

DHHS – Department of Health and Human Services

(3)

Filtering label

X-Ray Filters

fix: XX: X,X mm

XX: X,X mm

XX: X,X mm

switchable: XX: X,X mm

Refer to the original label that is applied on the device for information about the X-ray filters. This label is an example only.

(4)

Product identification label

SIEMENS		
		
(1P) Model No.	XXXXXXX	
		
(S) Serial No.	XXXX	IVK
Siemens AG, Wittelsbacherplatz 2, DE-80333 Muenchen		
Made in Germany		

– or –

Product identification label


SIEMENS		
XXXXXXXXXX		
	(01) XXXXXXXXXXXXXXXX GTIN	
	(240) XXXXXXXX MODEL	
	(422) XXX ORIGIN YY	IVK
	(21) XXXXX SERIAL	
	(20) XX REV	
Siemens Healthcare GmbH, Henkestr. 127, 91052 Erlangen, Germany		
Made in Germany		

(5)

Product status label

Rev.	00	01	02	03	04	05	06	07	08	09	10	11
12	13	14	15	16	17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32	33	34	35	36	37

High voltage generator

(1) Product status label

(2) Approval label (USA, DHHS)

(3) Product identification label

(4) Spare part identification label

(1)

Product status label

Rev.	00	01	02	03	04	05	06	07	08	09	10	11
12	13	14	15	16	17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32	33	34	35	36	37

Refer to the original label that is applied on the device for information about the current product version. This label is an example only.

(2)

Approval label (USA, DHHS)

This product complies with DHHS regulations
21 CFR Subchapter J, applicable at date of
manufacture.


Manufactured: XXXXXXXXXXXX

Siemens Healthcare GmbH
Henkestr. 127, 91052 Erlangen
Germany

DHHS – Department of Health and Human Services


(3)

Product identification label


(4)

Spare part identification label

Detector DMS A

(1) Product identification label

(1)

Product identification label


SIEMENS		
		
(1P) Model No.	XXXXXXX	
		
(S) Serial No.	XXXX	IVK
Siemens AG, Wittelsbacherplatz 2, DE-80333 Muenchen		
Made in Germany		

– or –

Product identification label

SIEMENS		
XXXXXXXX		
	(01) xxxxxxxxxxxxxxxx	GTIN
	(240) xxxxxxxx	MODEL
	(422) xxx	ORIGIN YY
	(21) xxxxx	SERIAL
	(20) xx	REV
IVK		
 Siemens Healthcare GmbH, Henkestr. 127, 91052 Erlangen, Germany		
Made in Germany		

Detector DMS B

(1) Product identification label

(1)

Product identification label

SIEMENS		
		
(1P) Model No.	XXXXXXX	
		
(S) Serial No.	XXXX	IVK
Siemens AG, Wittelsbacherplatz 2, DE-80333 Muenchen		
Made in Germany		


– or –

Product identification label

SIEMENS		
XXXXXXXX		
	(01) xxxxxxxxxxxxxxxxxxxx GTIN (240) xxxxxxxx MODEL (422) xxx ORIGIN YY (21) xxxxx SERIAL (20) xx REV	
IVK Siemens Healthcare GmbH, Henkestr. 127, 91052 Erlangen, Germany Made in Germany		

Patient table

Patient table with 2000 mm scan range

(1) Warning, crushing hazard for hand

(2) Caution, leg injuries

- (3) Approval label (USA and Canada, NRTL)
- (4) Duty cycle label
- (5) Approval label (USA, DHHS)
- (6) Product identification label (Patient table)
- (7) Type B applied part
- (8) Product identification label (Table top)
- (9) Total working load label
- (10) Equipotentiality label

(1)


Warning, crushing hazard for hand

Increased risk that hands can get pinched between table top and patient table.

(2)

Caution, leg injuries

Contusion of legs is possible when lowering the patient table.

Caution: Consult accompanying documents.

(3)

Approval label (USA and Canada, NRTL)

NRTL – Nationally Recognized Testing Laboratory

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(5)

Approval label (USA, DHHS)

This product complies with DHHS regulations
21 CFR Subchapter J, applicable at date of
manufacture.

Manufactured: XXXXXXXXXXXX

Siemens Healthcare GmbH
Henkestr. 127, 91052 Erlangen
Germany

DHHS – Department of Health and Human Services

(6)

Product identification label

SIEMENS			
XXXXXXXX			
	(01) XXXXXXXXXXXXXXXX	GTIN	
	(240) XXXXXXXX	MODEL	
	(422) XXX	ORIGIN YY	IVK
	(21) XXXXX	SERIAL	
 Siemens Healthcare GmbH, Henkestr. 127, 91052 Erlangen, Germany			
Made in Germany			

(7)

Type B applied part


Equipment is not suitable for direct cardiac application.

Indicates the degree of patient protection of this applied part (Type B applied part according to IEC 60601-1).

(8)

Product identification label

(9)

Total working load label

The maximum permissible working load must not exceed 227 kg (500 lbs), which is approximately 2226 N.


(10)

Equipotentiality label

Indicates the location of a pin that can be used for equalization with other medical devices in accordance with IEC 60601-1 and DIN 42801-2.

Multi-purpose table

- (1) Warning, crushing hazard for hand
- (2) Product identification label (Table top)
- (3) Equipotentiality label
- (4) Total working load label
- (5) Type B applied part
- (6) Approval label (USA and Canada, NRTL)
- (7) Duty cycle label
- (8) Product identification label (Multi-purpose table)
- (9) Approval label (USA, DHHS)
- (10) Caution, leg injuries


- (1)
Warning, crushing hazard for hand

Increased risk that hands can get pinched between table top and patient table.

(2)

Product identification label

(3)


Equipotentiality label

Indicates the location of a pin that can be used for equalization with other medical devices in accordance with IEC 60601-1 and DIN 42801-2.

(4)

Total working load label

The maximum permissible working load must not exceed 227 kg (500 lbs), which is approximately 2226 N.

– or –

Total working load label

The maximum permissible working load must not exceed 307 kg (677 lbs), which is approximately 3010 N.

(5)

Type B applied part


Equipment is not suitable for direct cardiac application.

Indicates the degree of patient protection of this applied part (Type B applied part according to IEC 60601-1).

(6)

Approval label (USA and Canada, NRTL)

NRTL – Nationally Recognized Testing Laboratory

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

(9)


Approval label (USA, DHHS)

DHHS – Department of Health and Human Services

(10)


Caution, leg injuries

Contusion of legs is possible when lowering the patient table.

Caution: Consult accompanying documents.

Physiological measurement module


- (1) Warning label
- (2) ECG socket
- (3) Type BF applied part
- (4) ECG channel selection

(1) Warning Label 	This equipment may be dangerous to patients and users unless safe expose factors, Instructions for Use, and maintenance schedules are observed. Indicates that the Type BF applied part is not defibrillation proof.
(2) ECG socket 	
(3) Type BF applied part 	Equipment is not suitable for direct cardiac application. Indicates the degree of patient protection of this applied part (Type BF applied part according to IEC 60601-1).
(4) ECG channel selection 	

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
- (1) Warning label
- (2) ECG socket
- (3) Type BF applied part
- (4) ECG channel selection

(1) Warning Label 	This equipment may be dangerous to patients and users unless safe expose factors, Instructions for Use, and maintenance schedules are observed. Indicates that the Type BF applied part is not defibrillation proof.
	(2) ECG socket 
(3) Type BF applied part 	Equipment is not suitable for direct cardiac application. Indicates the degree of patient protection of this applied part (Type BF applied part according to IEC 60601-1).
	(4) ECG channel selection 

Power distribution cabinet

- (1) Approval label (Saudi Arabia only, MOH)
- (2) Product status label
- (3) Product identification label
- (4) Input power label
- (5) Warning, check for absence of power

(1)

Approval label (Saudi Arabia only, MOH)

MOH – Ministry of Health

(2)

Product status label

Rev.	00	01	02	03	04	05	06	07	08	09	10	11
12	13	14	15	16	17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32	33	34	35	36	37

Refer to the original label that is applied on the device for information

about the current product version. This label is an example only.

(3)

Product identification label

XXXXXXXXXXXX

(01) xxxxxxxxxxxxxxxx

(240) xxxxxxxx

(422) xxx

(21) xxx

GTIN

MODEL

ORIGIN YY

SERIAL

IVK

(4)

Input power label

3N/PE~

380V/400V/420V/440V/460V/480V

50/60Hz

Continuous operation with intermittent loading

momentary operation / loading: XXX-XXxA

long time operation / loading: XX-XXA

(5)

Warning, check for absence of power

Warning

Even when the CT system is switched off there is voltage present. To ensure safe working conditions, all external live lines must be switched off on-site and secured against being turned on again. Before beginning work wait 2 minutes and check for absence of power.

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
6.5.2 Labeling Information and requirements VB20 version

Gantry

Front Side

(1) Approval (Saudi Arabia only, MOH)


MOH – Ministry of Health

(2) Follow the Instructions for Use

The user must consult and follow the Instructions for Use before using the CT system and its equipment. An electronic copy of the Instructions for Use is available at <http://www.siemens.com/medImaging-Manuals>

Rear side

- (1) Warning, check for absence of power
- (2) Product identification label (System)
- (3) Approval label (Malaysia only, MDA)
- (4) Product identification label (Gantry)
- (5) Approval label (USA only, DHHS)
- (6) Approval label (USA and Canada, NRTL)
- (7) Radio frequency information (Only for systems with an optional wireless i-Control)

(1) Warning, check for absence of power


Warning

Even when the CT system is switched off there is voltage present. To ensure safe working conditions, all external live lines must be switched off on-site and secured against being turned on again. Before beginning work wait 2 minutes and check for absence of power.

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(3) Approval label (Malaysia, MDA)

Registration no.: XXXXXXXXXXXX

Malaysia Authorised Representative:

Siemens Healthcare Sdn Bhd
Level 16-1 & 16-2 CP Tower,
No. 11, Jalan 16/11, Pusat Dagang Seksyen 16,
46350 Petaling Jaya
Phone: 1800 888 872

MDA – Medical Device Authority

(4) Product identification label

SIEMENS

XXXXXXXX

 (01) xxxxxxxxxxxxxxxx GTIN
(240) xxxxxxxx MODEL
(422) xxx ORIGIN YY IVK
(21) xxxxx SERIAL
(20) xx REV

Siemens Healthcare GmbH,
Henkestr. 127, 91052 Erlangen, Germany
Made in Germany

(5) Approval label (USA only, DHHS)

This product complies with DHHS regulations
21 CFR Subchapter J, applicable at date of
manufacture.
Manufactured: XXXXXXXXXXXX
Siemens Healthcare GmbH
Henkestr. 127, 91052 Erlangen
Germany


DHHS – Department of Health and Human Services

(6) Approval label (USA and Canada, NRTL)

NRTL – Nationally Recognized Testing Laboratory


(7) Radio frequency information

Only for systems with an optional wireless i-Control.

The individual, country specific approval may vary on this label.

Gantry frame

(1) Product status label

(2) Product identification label (Gantry)

(3) Approval label (USA only, DHHS)

(4) Approval label (USA and Canada, NRTL)


(5) Product identification label (System)

(1) Product status label

Rev.	00	01	02	03	04	05	06	07	08	09	10	11
	12	13	14	15	16	17	18	19	20	21	22	23
	24	25	26	27	28	29	30	31	32	33	34	35
	36	37										

Refer to the original label that is applied on the device for information about the current product version. This label is an example only.

(2) Product identification label

(3) Approval label (USA only, DHHS)


This product complies with DHHS regulations
21 CFR Subchapter J, applicable at date of
manufacture.
Manufactured: XXXXXXXXXXXX
Siemens Healthcare GmbH
Henkestr. 127, 91052 Erlangen
Germany

(4) Approval label (USA and Canada, NRTL)


NRTL – Nationally Recognized Testing Laboratory

(5) Product identification label (System)

Tube housing


- (1) Product status label
- (2) Focal spot size
- (3) Product identification label (Tube assembly)
- (4) Approval label (USA, DHHS)
- (5) Approval label (USA and Canada, NRTL)
- (6) Product identification label (X-ray tube)

(1) Product status label

Rev.	00	01	02	03	04	05	06	07	08	09	10	11
12	13	14	15	16	17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32	33	34	35	36	37

Refer to the original label that is applied on the device for information about the current product version. This label is an example only.


(2) Focal spot size

(3) Product identification label

(4) Approval label (USA, DHHS)

DHHS – Department of Health and Human Services

This product complies with DHHS regulations
21 CFR Subchapter J, applicable at date of
manufacture.

Manufactured: XXXXXXXXXXXX
Siemens Healthcare GmbH
Henkestr. 127, 91052 Erlangen
Germany

DHHS – Department of Health and Human Services

(5) Approval label (USA and Canada, NRTL)


NRTL – Nationally Recognized Testing Laboratory

(6) Product identification label

Tube collimator

- (1) RoHS label
- (2) Approval label (USA, DHHS)
- (3) Filtering label
- (4) Product identification label
- (5) Product status label


(1) RoHS label

(2) Approval label (USA, DHHS)


DHHS – Department of Health and Human Services

DHHS – Department of Health and Human Services

(3) Filtering label

Refer to the original label that is applied on the device for information about the X-ray filters. This label is an example only.

(4) Product identification label


(5) Product status label

Rev.	00	01	02	03	04	05	06	07	08	09	10	11
12	13	14	15	16	17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32	33	34	35	36	37

Refer to the original label that is applied on the device for information about the current product version. This label is an example only.

High voltage generator


- (1) Product status label
- (2) Approval label (USA, DHHS)
- (3) Product identification label
- (4) Spare part identification label

(1) Product status label

Rev.	00	01	02	03	04	05	06	07	08	09	10	11
12	13	14	15	16	17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32	33	34	35	36	37


Refer to the original label that is applied on the device for information about current product version. This label is an example only

(2) Approval label (USA, DHHS)


DHHS – Department of Health and Human Services


(3) Product identification label

(4) Spare part identification label


Detector DMS A


(1) Product identification label

(1) Product identification label


Detector DMS B


(1) Product identification label

- (1) Product identification label


-or-

Patient table

Patient table with 2000 mm scan range


- (1) Warning, crushing hazard for hand
- (2) Caution, leg injuries
- (3) Approval label (USA and Canada, NRTL)
- (4) Duty cycle label
- (5) Approval label (USA, DHHS)
- (6) Product identification label (Patient table)
- (7) Type B applied part
- (8) Product identification label (Table top)
- (9) Maximum load capacity label
- (10) Equipotentiality label

- (1) Warning, crushing hazard for hand

Increased risk that hands can get pinched between table top and patient table

- (2) Caution, leg injuries

Increased risk that body parts can get pinched between table and floor.

Take care when lowering the patient table to avoid injury to body parts.

To avoid damaging the patient table, make sure no accessories (for example, a trolley or cart) are placed underneath the patient table while lowering it.

- (3) Approval label (USA and Canada, NRTL)


NRTL – Nationally Recognized Testing Laboratory

- (4) Duty cycle label


The cyclic factor (c.d.f) specifies the ratio of the maximum operation time of the component to the total operation time of the system

- (5) Approval label (USA, DHHS)

DHHS – Department of Health and Human Services

(6) Product identification label


(7) Type B applied part

Equipment is not suitable for direct cardiac application. Indicates the degree patient protection of this applied part (Type B applied part according to IEC 60601-1).

(8) Product identification label

(9) Total working load label


The maximum permissible working load must not exceed 227 kg (500 lbs), which is approximately 2226 N.

(10) Equipotentiality label

Indicates the location of a pin that can be used for equalization with other medical devices in accordance with IEC 60601-1 and DIN 42801-2.

Multi purpose table


- (1) Warning, crushing hazard for hand
- (2) Product identification label (Table top)
- (3) Equipotentiality label
- (4) Maximum load capacity label
- (5) Type B applied part
- (6) Approval label (USA and Canada, NRTL)
- (7) Duty cycle label
- (8) Product identification label (Multi-purpose table)
- (9) Approval label (USA, DHHS)
- (10) Caution, leg injuries

- (1) Warning, crushing hazard for hand

Increased risk that hands can get pinched between table top and patient table

- (2) Product identification label


(3) Equipotentiality label

Indicates the location of a pin that can be used for equalization with other medical devices in accordance with IEC 60601-1 and DIN 42801-2.

(4) Maximum load capacity

The maximum permissible working load must not exceed 227 kg (500 lbs), which is approximately 2226 N.

-or-

The maximum permissible working load must not exceed 307 kg (667 lbs), which is approximately 3010 N.

(5) Type B applied part

Equipment is not suitable for direct cardiac application.

Indicates the degree of patient protection of this applied part (Type B applied part according to IEC 60601-1).

(6) Approval label (USA and Canada, NRTL)


NRTL – Nationally Recognized Testing Laboratory

(7) Duty cycle label


The cyclic duration factor (c.d.f) specifies the ratio of the maximum operation time of the components to the total operation time of the system.

(8) Product identification label


(9) Approval label (USA, DHHS)

DHHS – Department of Health and Human Services

(10) Caution, leg injuries


Increased risk that body parts can get pinched between table and floor.

Take care when lowering the patient table to avoid injury to body parts.

To avoid damaging the patient table, make sure no accessories (for example, a trolley or cart) are placed underneath the patient table while lowering it.

Physiological measurement module

- (1) Warning label
- (2) ECG socket
- (3) Type BF applied part
- (4) ECG channel selection

- (1) Warning label

This equipment may be dangerous to patients and user unless safe expose factors, Instructions for Use, and maintenance schedules are observed.

The Type BF applied part is not defibrillation proof.

- (2) ECG socket

- (3) Type BF applied part


Equipment is not suitable for direct cardiac application.

Indicates the degree of patient protection of this applied part (Type BF applied part according to IEC 60601-1).

- (4) ECG channel selection

Physiological measurement module for multipurpose table

- (1) Warning label
- (2) ECG socket
- (3) Type BF applied part
- (4) ECG channel selection

(1) Warning label

This equipment may be dangerous to patient and user unless safe expose factors, Instructions for Use, and maintenance schedules are observed.

Indicates that the Type BF applied part is not defibrillation proof.

(2) ECG socket

(3) Type BF applied part

Equipment is not suitable for direct cardiac application.

Indicates the degree of patient protection of this applied part (Type BF applied part according to IEC 60601-1).


(4) ECG channel selection

6.5.3 Labeling Information and requirements VA48 version

Gantry

Front side

1) Follow Instruction for Use

Follow Instruction for Use


Obliges the user to consult the Instructions for Use before using the CT system.

An electronic copy of the Instructions for Use is available at <http://www.siemens.com/MedImaging-Manuals>.

(2) Ministry of Health, Kingdom of Saudi Arabia
(Saudi Arabia only)

Rear side

SIEMENS	
SOMATOM XXXXXXXX	
(1P) Model No.	XXXXXXX
(S) Serial No.	XXXX
Siemens Aktiengesellschaft, Wittelsbacher Platz 2, DE-80333 Muenchen	
Made in Germany	

(1) Manufacturer identification label (System)

SIEMENS	
XXXXXXX	
(1P) Model No.	XXXXXXX
(S) Serial No.	XXXX
Siemens AG, Wittelsbacherplatz 2, DE-80333 Muenchen	
Made in Germany	

(2) Manufacturer identification label (Gantry)

This product complies with DHHS regulations
21 CFR Subchapter J, applicable at date of
manufacture.

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

(1) Manufacturer identification label (Gantry)

This product complies with DHHS regulations
21 CFR Subchapter J, applicable at date of
manufacture.

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

(4) Manufacturer identification label (System)

Rev.	00	01	02	03	04	05	06	07	08	09	10	11
12	13	14	15	16	17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32	33	34	35	36	37

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This product complies with DHHS regulations
21 CFR Subcpter J, applicable at date of
manufacture.

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
(3) Manufacturer identification label

(4) Focal spot size



(5) Nationally Recognized Testing Laboratory (an organization which is recognized by OSHA e.g.: UL, CSA, TÜV Süd)

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
(2) Filtering label

(3)

This product complies with DHHS regulations
21 CFR Subchapter J, applicable at date of
manufacture.

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


Mat. No. XXX XX XXX
Rev 00 01 02 03 04 05 06

(3) Product status label

**This product complies with DHHS regulations
21 CFR Subcpter J, applicable at date of
manufacture.**

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

(2) Label for X-ray products


Department of Health and Human Services (DHHS)

(3) Nationally Recognized Testing Laboratory (an organization which is recognized by OSHA e.g.: UL, CSA, TÜV Süd)

– or –


(4) Total working load label

The label is depicted on the positioning aids to indicate the maximum loads.

(5) Degree of protection label type B

Type B applied part (degree of protection against electrical shock in accordance with IEC 60601-1).

(6) Warning label of leg injuries

Contusion of legs is possible when moving the patient table.

(7) Duty cycle label

The cyclic duration factor (c.d.f) specifies the ratio of the maximum operation time of the component to the total operation time of the system.

(8) Equipotential bonding label

The label indicates the location of a pin for potential equalization with other medical devices in accordance with IEC 60601-1 and DIN 42801-2.

(9) Warning label of electric shock

Risk of electric shock



Refer servicing to qualified service personnel. To avoid risk of injury from electric shock, do not open this enclosure.

(10) Warning label of hand injuries

Increased risk that hands can get pinched between flat table top and patient table.

Physiological measurement module

(1) Operating instructions

Symbol explanation:

Label: ECG channel selection



Label: ECG socket

Label: Type BF applied part (in accordance with IEC 60601-1)

Physiological measurement module for multi-purpose table

(1) Operating instructions

Symbol explanation:

Label: ECG channel selection


Label: ECG socket


Label: Type BF applied part (in accordance with IEC 60601-1)

Power distribution cabinet

(1) Manufacturer identification label

(2) Warning (check de-energized)

(3) Ministry of Health, Kingdom of Saudi Arabia
(Saudi Arabia only)

(4) Product status label

6.6 Artwork of the MD labeling in Russian

Additionally, a Russian label is applied with the following content (please note that Labeling mentioned below is applied only for the SOMATOM Definition Flash systems which are delivered to Russia not earlier than September 2020):

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7. Sterility (if applicable)

The medical device is non-sterile.

8. The list of International Standards/ Directives/ Documents applicable to the Medical Device

Certificates

The CT system has the following certificates:

- EC Declaration of Conformity according to Medical Device Directive 93/42/EEC
- NRTL/CSA-Certificate according to ANSI/AAMI ES60601-1:2005/

(R)2012 and A1:2012 and CAN/CSA-C22.2 NO. 60601-1:2014

Accessory equipment connected to the analog and digital interfaces must be certified according to the respective IEC standards. Furthermore, all configurations must comply with the applicable standards, for example, IEC 60950 for data processing equipment and IEC 60601-1 for medical electrical equipment.

Any person who connects additional equipment to the signal input part or signal output part is configuring a medical system and is therefore responsible for ensuring that the system complies with the applicable standards.

Statutory regulations for medical device

The CT system and its components comply with the following major statutory regulations:

Statutory regulation	Title
93/42/EEC	Council Directive 93/42/EEC of 14 June 1993 concerning medical devices
2014/53/EU	DIRECTIVE OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 April 2014 on the harmonization of the laws of the Member States relating to the making available on the market of radio equipment and repealing Directive 1999/5/EC
21 CFR 801.109(b)(1)	Prescription devices
21 CFR 1020.30	Diagnostic x-ray systems and their major components
21 CFR 1020.33	Computed tomography (CT) equipment
21 CFR 1040.10	Laser products
21 CFR 1040.11	Specific purpose laser products

Product environmental protection regulations

The CT system and its components comply with the following product environmental protection regulations.

Standard	Title
2006/66/EC	Directive 2006/66/EC of the European Parliament and of the Council of 6 September 2006 on batteries and accumulators and waste batteries and accumulators and repealing Directive 91/157/EEC
2011/65/EU	Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment
1907/2006/EC	Regulation (EC) No 1907/2006 of the European Parliament and of the Council of 18 December 2006 concerning the Registration, Evaluation, Authorization and Restriction of Chemicals (REACH), establishing a European Chemicals Agency, amending Directive 1999/45/EC and repealing Council Regulation (EEC) No 793/93

and Commission Regulation (EC) No 1488/94 as well as Council Directive 76/769/EEC and Commission Directives 91/155/EEC, 93/67/EEC, 93/105/EC and 2000/21/EC

Standards according to IEC 60601-1, 3rd edition

Standard	Title
IEC 60601-1:2005+A1:2012	Medical electrical equipment – Part 1: General requirements for basic safety and essential performance
IEC 60601-1-2:2007	Medical electrical equipment –
IEC 60601-1-2:2014	Part 1-2: General requirements for basic safety and essential performance –Collateral standard: Electromagnetic compatibility – Requirements and tests
IEC 60601-1-3:2008+A1:2013	Medical electrical equipment – Part 1-3 : General requirements for basic safety and essential performance –Collateral standard: Radiation protection in diagnostic X-ray equipment
IEC 60601-1-6:2006	Medical electrical equipment –
IEC 60601-1-6:2010+A1:2013	Part 1-6: General requirements for basic safety and essential performance –Collateral standard: Usability
IEC 60601-2-28:2010	Medical electrical equipment – Part 2-28: Particular requirements for the basic safety and essential performance of X-ray tube assemblies for medical diagnosis
IEC 60601-2-44:2009+A1:2012	Medical electrical equipment – Part 2-44: Particular requirements for the basic safety and essential performance of X-ray equipment for computed tomography

Standards in the scope of IEC 60601-1

Standard	Title
IEC 60825-1:2007	Safety of laser products –
IEC 60825-1:2014	Part 1: Equipment classification and requirements
IEC 61223-2-6:2006	Evaluation and routine testing in medical imaging departments –
	Part 2-6: Constancy tests – Imaging performance of computed tomography X-ray equipment
IEC 61223-3-5:2004	Evaluation and routine testing in medical imaging departments –Part 3-5: Acceptance tests –
	Imaging performance of computed tomography X-ray equipment
IEC 62304:2006+A1:2015	Medical device software – Software life cycle processes
IEC 62366:2007	Medical devices - Application of usability engineering to medical devices
IEC 62366:2007+A1:2014	
EN 50419:2005	Marking of electrical and electronic equipment in accordance with Article 11(2) of Directive 2002/96/EC (WEEE)
ISO 10993-1:2009	Biological evaluation of medical devices –
	Part 1: Evaluation and testing within a risk management process

Referenced standards

The instructions for use refer to the following standards:

Standard	Title
DIN 42801-2:1984	Connection device for potential equalization conductors
DIN 5032-7:2017	Photometry – Part 7: Classification of illuminance meters and luminance meters
DIN 6868-4:2007	Image quality assurance in diagnostic X-ray departments – Part 4: Constancy testing of medical X-ray equipment for fluoroscopy

DIN 6868-157:2014	Image quality assurance in diagnostic X-ray departments – Part 157: X-ray ordinance acceptance and constancy test of image display systems in their environment
EN-50419:2006	Marking of electrical and electronic equipment in accordance with Article 11(2) of Directive 2002/96/EC (WEEE)
IEC-60336:2005	Medical electrical equipment – X-ray tube assemblies for medical diagnosis – Characteristics of focal spots
IEC-60417-DB-HS:2008	Graphical symbols for use on medical equipment – Part 1: General
IEC 60529:1989+A1:1999+A2:2013	Classification of degrees of protection provided by enclosures
IEC 60601-2-7:1998	Medical electrical equipment – Part 2-7: Particular requirements for the safety of high-voltage generators of diagnostic X-ray generators
IEC 60721-3-2:1997	Classification of environmental conditions – Part 3: Classification of groups of environmental parameters and their severities – Section 2: Transportation
IEC 60950:2013	Information technology equipment – Safety – Part1: General requirements
IEC 61000-3-2:2009	Electromagnetic compatibility (EMC) – Part 3-2: Limits - Limits for harmonic current emissions (equipment input current ≤16 A per phase)
IEC 61000-3-3:2013	Electromagnetic compatibility (EMC) – Part 3-5: Limits - Limitation of voltage fluctuations and flicker in low-voltage power supply systems for equipment with rated current greater than 75 A
IEC 61000-4-2:2008	Electromagnetic compatibility (EMC) – Part 4-2: Testing and measuring techniques – Electrostatic discharge immunity test

IEC 61000-4-3:2010	Electromagnetic compatibility (EMC) – Part 4-3: Testing and measurement techniques – Radiated, radio-frequency, electromagnetic field immunity test
IEC 61000-4-4:2012	Electromagnetic compatibility (EMC) – Part 4-4: Testing and measuring techniques – Electrical fast transient/burst immunity test
IEC 61000-4-5:2005	Electromagnetic compatibility (EMC) – Part 4-5: Testing and measurement techniques – Surge immunity test
IEC 61000-4-6:2013	Electromagnetic compatibility (EMC) – Part 4-6: Testing and measurement techniques – Immunity to conducted disturbances, induced by radio-frequency fields
IEC 61000-4-8:2009	Electromagnetic compatibility (EMC) – Part 4-8: Testing and measurement techniques – Power frequency magnetic field immunity test
IEC 61000-4-11:2004	Electromagnetic compatibility (EMC) – Part 4-11: Testing and measurement techniques – Voltage dips, short interruptions and voltage variations immunity tests
IEC 61223-2-5	Evaluation and routine testing in medical imaging departments – Part 2-5: Constancy tests, Image display devices
IEC 61267:2005	Radiation conditions for use in the determination of characteristics
IEC 61674:2012	Medical electrical equipment – Dosimeters with ionization chambers and/or semiconductor detectors as used in X-ray diagnostic imaging.
ISO 4037-4:2004	X and gamma reference radiation for calibrating dosimeters and doserate meters and for determining their response as a function of photon energy.

Standards referenced in the tube manual

The listed standards are referenced in the Straton MX P Instructions for Use.

Standard	Title
EN 50581:2012	Technical documentation for the assessment and electronic products with respect to the restriction of hazardous substances.
IEC 60336:2005+C1:2006	Medical electrical equipment – X-ray tube assemblies for medical diagnosis – Characteristics of focal spots
IEC 60336:2006	Medical electrical equipment – X-ray tube assemblies for medical diagnosis – Characteristics of focal spots
IEC 60522:1999	Determination of the permanent filtration of X-ray tube assemblies
IEC 60522:2003	Determination of the permanent filtration of X-ray tube assemblies
IEC 60526:1978	High-voltage cable plug and socket connections for medical X-ray equipment
IEC 60526:1978+C1+2010	High-voltage cable plug and socket connections for medical X-ray equipment
IEC 60601-1:2005+C1:2006+C2:2007	Medical electrical equipment – Part 1: General requirements for basic safety and essential performance
IEC 60613:1989	Electric, thermal and loading characteristics of rotating anode X-ray tubes for medical diagnosis
IEC 60613:2010	Electrical and loading characteristics of X-ray tube assemblies for medical diagnosis

Standards addressed to the system owner

The following standards describe test procedures and the parameters to be tested by the system owner. The CT system provides tools (software) to perform and evaluate the measurements required by the standards.

Standard	Title
DIN 6868-157:2014	Image quality assurance in diagnostic X-ray departments – Part 157: X-ray ordinance and constancy test of image display systems in their environment.
IEC 61223-2-4:1994	Evaluation and routine in medical imaging departments – Part 2-4: Constancy tests – Hard copy cameras
IEC 8001-1:2010	Application of risk management for IT-networks incorporating medical devices – Part 1: Roles, responsibilities and activities

9. List of equipment and measurement tools necessary for quality control of the medical device (if applicable)

Accessories for quality measurement

You require the set of phantoms for the quality measurement and possibly additional accessories for the constancy test.

Accessories for daily quality measurement

The following components are required for the daily quality measurement:

- Phantom holder
- Set of phantoms

Accessories for the monthly constancy test

For the monthly constancy test, you require the following phantoms and accessories:

- Phantom holder
- Quality phantom
- 100 kg (220 lbs) weight (a patient equivalent load not exceeding 135 kg (297 lbs))
- Ruler, 40 cm or longer

Accessories for the Low Contrast test

For the low contrast test, you require the following phantoms and accessories:

- Phantom holder
- Low contrast phantom

10. Transportation conditions and requirements


10.1. Acceleration

To avoid system damage, the following conditions are mandatory for transport

- Maximum acceleration according to DIN EN 60721-3-2, Class 2M2.

Directions of impact:

Fig. 30: Directions of impact

(X) Permissible impact: 4 g (corresponds to shunting a train)

(Z) Permissible impact: 2 g

Transport all CT components in the respective transport packaging provided by the manufacturer. When transporting the CT system across country terrain, only use trucks with air shock suspension. **Additionally, observe the "Safety information" section in the first chapter "General information".**

10.2. Climatic conditions for transport and storage

10.2.1. Standard climatic conditions

The climatic requirements for storing the CT system are identical to the climatic requirements for transportation.

The data for transporting and storing the CT system are applicable only if the system is shipped free of damage in the transport packaging provided by the manufacturer.

Tab. 58 Climatic conditions: Transport and storage


Designation	Range
Temperature	-20 °C to +50 °C
Relative humidity	10% to 90%
Temperature gradient	10 K/h
Barometric pressure	700 hPa to 1060 hPa
Max. storage time	2 months

All components have to be stored in the respective transport packaging provided by the manufacturer.

Use equivalent packaging for intermediate storage if the original packaging provided by the manufacturer is no longer available.

Fig. 72: Transport and storage conditions (climatogram)

After moving the CT system from a cold environment into a warm room, make sure that the entire CT system is dry before switching it ON.

10.2.2. Extended climatic conditions

If the temperature range during transport and storage exceeds the specified values

- from -20 °C to -45 °C and
- from 50 °C to 70 °C

special packaging must be used; it limits the temperature for the system to the specified values. Verify that the components to be transported can withstand these temperatures without damage.

11.Storage and Shelf Life

Climatic conditions for storage:

Designation	Range
Temperature	-20 °C to +50 °C
Relative humidity	10% to 90%
Temperature gradient	10 K/h
Barometric pressure	700 hPa to 1060 hPa
Max. storage time	2 months

Service life: 10 years

12.Procedure of recycling and disposal (Medical Device recycling or disposal data)

Current recycling and waste management regulations assign product responsibility to the manufacturer / seller. Compliance with these product regulations includes evaluation of recyclables and environmentally safe disposal of materials and components.

This document provides information about the type, quantity and location of hazardous materials in the as well as instructions how to remove the relevant components. Follow all guidelines given in the Disposal Instructions when disposing of the system. When disposing of the system and/or parts, please note that repairable parts have to be returned for reuse.

CAUTION

[hz_serdoc_F14G08U01M02]

Hazardous materials have to be properly disposed of by an authorized waste management company.

Risk of environmental pollution!

⇒ **Dispose of hazardous materials in accordance with local regulations.**

13.Information about maintenance (general and on the territory of Russia). Requirements to the authorized personnel who it carries. List of recommended tools and instruments for maintenance and repair. Requirements to the maintenance intervals (if necessary):

On request, you can have maintenance work performed regularly by the Siemens Customer Service. Please contact your Siemens Service to agree on a maintenance contract for your system.

Service tools

SOMATOM Definition Flash

Auxiliary materials – Definition

05729496

Part-No.	Description
----------	-------------

Comment

03473287	GREASE TUNGREASEBS =G=
----------	------------------------

For lubrication of the lifting s pindle on the Z4-treatment couch.

04007725	GREASE TOPASNCA52 =G=
----------	-----------------------

03499894	GREASE MOBILTEMPSHC100 – 380g =G=
----------	-----------------------------------

Mobiltemp SHC 100 is recommended for use in high speed bearings where a wide temperature range is desired, and also in high temperature applications where extended service between relubrication intervals is required. It is particularly suitable for use in electric motor bearings where operating conditions demand reduced friction, low wear and long service life. Mobiltemp SHC 100 has recommended operating temperature range of 50°C to 200°C with appropriate relubrication intervals. 380g Cartridge.

02048866

SCREW SEALANT Loctite 270 =G=

LOCTITE 270 is a high-strength threadlocker for maximum efficacy in the securing and sealing of bolts, nuts and studs to prevent loosening due to vibration. The product serves to permanently lock assemblies which must not come loose.

02048874

SCREW SEALANT Loctite 221 =G=

LOCTITE 221 Low strength, low viscosity, small threads

Threadlocker, Breakaway torque 8.5Nm

05507038	SCREW SEALANT Loctite 243 =G=
----------	-------------------------------

Medium strength adhesive/sealant to lock and seal nuts, bolts, screws and fittings against vibration loosening leakage and thread corrosion. Part can be disassembled with hand tools.

04009648	PNT. Marker Hellgrün P0.05.75 =G=
----------	-----------------------------------

Light green

02884609	HIGH PRESS. GREASE GUN
----------	------------------------

Comes with the grease nipple type K/H and D, also one extension pipe 3,9" length.
For lubrication of:

1. The spindle on SOMATOM, filled with Unitemp 500 grease.
2. The ball bearings of the Polystar, filled with Mobiltemp SHC 100 (P/N: 34 99 894)
3. The SG3, filled with lithium soap grease Calypsol H 443
4. The magnetic brace at Siretom 1, filled with Optimoly.

Consequently, 4 grease guns are required. Contents: 150 ccm

10569439 Transport compensation container

For Part-No. 07740769, Definiton only!

10610161

Transport compensation container

Only for

Def. AS, Edge, Flash and SOMATOM Drive

10545450 PRIMER CMP 1K GRAU =G=

10545451 ADHES. Loctite E18 =G=

03804676 Large grease gun

07757524 Nozzle pipe with peak

Not valid for CT

Nozzle pipe peak

07757516 Manual grease press

Not valid for CT!

The manual grease press

(filled with Mobiltemp SHC 100)

Is used for Urospkop ACCESS rails

And carries:

e.g. tabletop longitudinal,

tabletop transverse, I.I. carrier,

x-ray tube carrier.

(Partnumber 7757516 is without nozzle pipe with peak)

SOMATOM Definition Flash**Electrical and mechanical parts****10496938**

Part- No.	Description Comment
(04299 850)	Tool case, basic & suppl., magnetic
099006 63	Torx key set, L-shaped, TX6 to TX25 7-Piece Torx key set. The L-shaped wrenches includes TX 7, 8, 9, 10, 15, 20 and 25 in a plastic holder
034245 61	TORQUE WRENCH 25-130NMX1/2" The torque wrench (20- 95 ft.lb) has a accuracy of $\pm 4\%$ from the preset value. No dials to watch-clicks when torque is reached. MAGNETIC.
07074594	TORQUE WRENCH 1-6NMX 1/4" The torque wrench (1-6 Nm) has an accuracy of $\pm 4\%$ from the present value
034300 63	TORQUE WRENCH 1-12NMX1/4" =D= For service and maintenance. The torque wrench has 2 scales, 2-12 Nm and 10-100 lb/in 1/4. Accuracy $\pm 4\%$ from the preset value. No dials to watch-clicks when torque is reachet.
(07059 975)	TORQUE WRENCH 2-20NMX1/4"
044151 13	SPRG BALAN 200 N For measuring tension- and traction forces. Scale division 4 N. Overload protection up to 20% over max. measuring value. Accuracy $\pm 1\%$ from max. scale reading. Calibration cycle: 2 years
080382 59	FEELER GAUGE 0,05MM-1MM For adjustment the Somatom collimator System. Set consists 18 blades from 0,05..0,10 to 2,00 mm.
097026 06	TKD FIELD SERVICE KIT This compact, portable kit quickly and reliable removes static charges on the service technician and also provides a static-free surface on wich to lay parts during service calls.

The kit includes two Charge-Guard wrist bands and one combination cord with snap on to the work surface. The longer cord connects the wrist band and the shorter cord attaches with an alligator clip to any reliable ground.

Each cord has its own independent 1 Mega Ohm resistor for user safety.

034684 GLOVE INDUSTR.HERREN –BAUMWOLLE

77

Gloves for servicing sliding-rings and MOS- PCB#s.

Size: One size fits all

101618 Balancing-set

72

Transportation device

077408 Transport frame_Matrix

43

for transportation of Gantry

113379 Transport frame DS cpl

28

Includes transport frame 10743385

101387

48

Transport frame_cpl 8/2

for transportation of PHS-4 and PHS-5

Local supply,

no Med Part-no. is available

7mm hex bit socket ,1/2" Service on Gantry

17mm hex bit socket, 1/2", long type Service on Gantry

1/2" extension, length >300mm Service on Gantry

3/8" extension, length 300mm Service on Gantry

2 water hoses with 1/2" connections Service on WCS

Fork wrench 41mm Connecting of water hoses

Fork wrench 46mm Service on WCS

Pipe wrench Service on WCS

Set of torx sockets 1/4" TX10-TX25

Set of torx sockets 3/8" TX30-TX50

Set of Allen sockets for torque wrench (size 3 - 14 metric) for Service

Set of hex bit sockets, 1/4" or 1/2" for torque wrench (size 3 - 17 metric) for Service

Socket adapter 1/4" drive - 1/2" fit for Service

Socket adapter 1/2" drive - 1/4" fit for Service

Socket adapter 1/2" drive - 3/8" fit

Cable ties Replacement of parts

Vacuum cleaner Maintenance

Service Laptop

for Service and documentation

Suction lifting device type 602.2 (floor tile remover)

Removal and Installation of Gantry covers

Tool case, basic suppl., magnetic

04299850

Part-No.	Description
----------	-------------

Comment

09702457	TOOL CASE, BASIC UNIT
----------	-----------------------

The toolcase contains 68 items and is almost similar to the former CT-tool case. This type of tool case is a basic kit and can be extended with the tool set for service engineers
Part No: 09702465 RE999.

09702465	SUPPLEMENTARY TOOL KIT
----------	------------------------

TORQUE WRENCH 5-25NMX1/4"

07059975

Part-No.	Description
----------	-------------

Comment

07059975	TORQUE WRENCH 5-25NMX1/4"
----------	---------------------------

The torque wrench has a accuracy
of $\pm 4\%$ from the preset value.
No dials to watch - clicks when
torque is reached.

SOMATOM Definition Flash

Measuring and testing tools

10188371

Part- No.	Description Comment
097027 13	ROTATING-FIELD INDICATOR For determine the phase sequence of three –phase AC-systems and to indicate the phase voltages with LED?s. According to DIN 57 413
102666 71	Endoscope camera Please only order this camera when performing UL CT076/14/P (Field inspection of monitor ceiling sys.)
099948 31	Fluke 287- Digital multimeter Calibration cycle: 1 year For all business units and systems. True rms Digital Multimeter includes:Hardware Users Manual 4 AA Batteries Silicone Test Leads (TL71). Replaces Fluke 8060A and Fluke 87. DC Voltage: 0 to 1000V DC Current: 0 to 10A (20A for 30 seconds) AC Voltage, True-rms: 2.5 mV to 1000V - 100 kHz bandwidth AC Current, True-rms: 25 mA to 10A (20A for 30 seconds) Resistance: 0 to 500 Mohms Conductance: 0 to 500 Nanosiemens Capacitance: 0.001 nF to 50 mF, Diode Test: 3.1V Temperature: -200°C to 1350°C (-328°F to 2462°F) Frequency: 0.5 Hz to 1000 kHz Accuracy(Basic DC V): 0.025% Accuracy(Basis AC V): 0.4%
107594 56	Fluke Scopemeter medical 190-202
(07752 848)	SPOTMETER fuer SMfit ACT =D=

070632 GONIOMETER NEIGUNG DIGITAL
90

073556 Belt,Tension Meter, Conti VSM-1
42

086142 Belt-tension meter
03

016217 HV DUMMY PLUGS E BEZ
91
Only used in HVT incl. 2pcs, Dummy plugs 3pol.

100937 Dummy Plugs
00
Only used in HV-sockets
Delivery unit contains 1 plug!

100959 Dosisfilm GAFCHROMIC XR-QA / XR-QA2
91
Total shelf life 365 days

034137 DIALGAUGE A1 DIN878-
88
Two pieces required for collimator adjustment on SOMATOM
Resolution 0,01 mm

-
Order directly:
Local supply;
No Med Part-no. is available

"Measurement device for protective ground measurement according to IEC62353,
respectively DIN VDE 0751.
The Secutest S3 shown in the image is only an example."

Measuring tape, metric
Adjusting the membrane position of the tube cooling system (Acron Q) and
Acceptance test

Dose meter Dali CT, Capintec WK92 or DIADOS
for Acceptance test

CCT champer Capintec PC4P with cable and holder
for Acceptance test

X-Ray film (13 x 18cm)

for Acceptance test

Spot-Luminance-Meter LXcan Please order this tool via: Weblogsx

Color Checker X-Rite M50101 Please order this tool via: Weblogsx

Inspection label for DIN 6868-157 Please order this label via: Weblogsx

SPOTMETER fuer SMfit ACT =D=

07752848

Part-No.	Description
	Comment
07752848	SPOTMETER fuer SMfit ACT =D= The device is used to Control/(Adjust CRT and TFT Monitors. Attention: The "old" MAVO Monitor may only be Used to Control/Adjust CRT Monitors
0864297	KABEL SMfit ACT bus cable kit

SOMATOM Definition Flash

10496937

Phantoms

Part-No.	Description
	Comment
10430785	Acceptance Set cpl., incl. Water phantom 300, cpl., with slice profile phantom
07032998	CLAMP LEVER 47MM LANG,SCHWARZ
07206006	PHANTOM HOLDER FOR SF6 Holder for the Paine phantom (USA only). Used on SOMATOM
09660663	Low-contrast phantom For measuring the low contrast resolution.

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Comment

09900788

Sticker / Complaint
Combined sticker/label to mark complaint components.

09900804

LABEL E (1 Pack = 100 pcs.)

For USA only!

To mark and seal SIEMENS-original parts. The label will be destroyed during removing
(New design!) One pack contains 100 labels.

09994773

RETURN GOODS INFO LABEL

1 Set contains 100 pcs.

Five-point, perforated label

"Return goods Information" including RD-note and error code. Basic color green.

For complete information please see reverse side of the label.

08428344

WARNING PLATE

For radiation protection/high voltage. use due to lokal regulation as required.

07281301

DEFECTIVE-LABEL FOR RETURN PARTS

10 Sheets of labels.

05551333

kit warning labels

Order directly:

Label for Monitor Acceptance test for Monitor

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Installation Template P47A

10430758

- Routine checks:
 - Image quality tests: Daily
 - Camera test: Daily
 - Function checks: Daily
 - Constancy test * : Monthly

*) We recommend that the constancy test be performed by Siemens service personnel in conjunction with a maintenance contract.
- Periodic maintenance activities are subdivided into:
 - Partial maintenance: Every 6 months
 - Full maintenance: Every 12 months
 - Preventive replacement of safety-relevant parts: Every 24 months

This work must be performed by qualified and authorized service engineers. In this context, qualified means that the service engineer has been trained and has had practical experience in the required routines so that he/she is capable of performing service work on the system. Authorized means that the service engineer is recognized by the system owner as a qualified service engineer and is therefore authorized to perform maintenance work on the system.

Upon system start-up, we recommend designating one employee who is responsible for ensuring that routine checks and preventive inspection and maintenance work are performed. The same employee should also be responsible for archiving all protocols in the "System Owner Manual".

In addition to repair services, Siemens Healthineers also offers its customers complete maintenance service, including inspection and preventive maintenance on Siemens Healthineers systems. These services can be provided either on a service call basis or in the form of a flexible maintenance contract. If you have not received a quotation for the performance of preventive maintenance and inspection from our Siemens Healthineers UPTIME Services, your local Siemens Healthineers sales representative will gladly provide this information for you.

On the territory of Russian Federation please use following contacts:

Siemens Healthcare LLC

Dubiniskaja str. 96

115093 Moscow, Russia

Tel.: +7 495 737 10 00

14. Information about calibration of the MD / requirements to calibration (if applicable)

System related artifacts

System related artifacts may be due to the following:

- Uncalibrated system
- Detector deviations

 **CAUTION**

When you restart the system, the detector has not yet reached operating temperature!

Wrong diagnosis due to image artifacts.

- ◆ Calibrate the system as part of the checkup. Repeat calibration if ring artifacts occur.

Calibration

If the CT system is not properly set up (for example, no checkup performed) CT scale displacements and inhomogeneities can result. This may result just after switching on the unit, during warm-up to operational temperature (calibration), or due to an extended length of service of the X-ray tube. The same applies to defects in the measurement system.

The daily quality checkup usually detects and reports such discrepancies. Other system defects or calibration deviations cause streaks or (partial) rings to appear in the CT image.

Detector deviations

In CT systems with a rotating combined tube-detector system, even minimal deviations of individual detector channels from the original calibration level may lead to rings or partial ring artifact structures in the CT image. The closer such channels lie to the detector center, the greater are these effects.

 **CAUTION**

When you restart the system, the detector has not yet reached operating temperature!

Wrong diagnosis due to image artifacts.

- ◆ Calibrate the system as part of the checkup. Repeat calibration if ring artifacts occur.

In an extreme case, rings may appear concentrated as a blurred spot in the center of the scan field (the scan field center appears at the image center if the Center X and Y image parameters are both zero). Such rings are easily recognizable, however, a blurred spot in the center of the scan field might lead to inaccurate diagnoses.

Repeating scans

Measurement related or system related artifacts cannot always be excluded. Therefore the scan should be repeated if a tomographic structure appears questionable. This should be done after slightly shifting the patient's position within the slice, e.g., by changing vertical position of the table by more than 5 mm. The result can be used to exclude system-dependent long-term or temporary equipment faults as the artifact cause.

Persisting object related artifacts can be considered measurement related with a high probability. They should be recognizable as such by an experienced radiologist.

15. Warranty

EXCERPT from the General Conditions of Supply and Delivery – Exports, Healthcare Divisions Imaging & Therapy Systems and Clinical Products engl_lieferbedingungen_Healthcare.dot (page 2)

10. Warranty

10.1. We warrant that the deliveries are free from material defects and lack none of the warranted characteristics at the time of risk transfer. Warranted characteristics are only those characteristics explicitly designated as such in the specifications. Our relevant Special Supply and Warranty Conditions shall apply in respect of products for which varying warranty regulations are agreed for product-specific reasons, such as high-vacuum elements or flat detectors.

10.2. The warranty period is 12 months from installation and commissioning, but shall expire no later than 15 months after delivery in accordance with the agreed Incoterms. The period of warranty for defects in respect of spare parts purchased from us is 6 months from the date of risk transfer.

10.3. Warranty claims shall not arise from trivial deviation from the agreed state or condition, in the case of non-reproducible software errors or in the course of normal wear and tear. Neither shall they arise from damages resulting from improper or negligent handling, assembly or commissioning by the Purchaser or third parties, excessive demand, improper use, failure to comply with operating instructions, use of inappropriate operating materials, incorrect structural prerequisites or technical specifications of the Purchaser or third parties, or due to special external influences not provided for in the contract. If modifications or repairs are performed by the Purchaser or a third party, no warranty claims shall exist in respect of such modifications or repairs or the consequences arising therefrom. Unless otherwise agreed, we accept no liability for material defects in used products.

10.4. Immediately after receiving the deliveries, the Purchaser shall examine them for defects. Obvious defects shall be reported immediately after shipment or installation or at the latest within fourteen (14) calendar days of receipt of the deliveries. Hidden defects shall be reported as soon as they are identified. Complaints shall be raised in writing, specifying the order data, the invoice and shipment numbers, and the identification number of the defective deliveries. If the Purchaser fails to provide notification in the correct form and in good time, the deliveries shall be deemed approved, subject to the reservation that hidden defects may be reported at a later time. The Purchaser shall not refuse to accept delivery due to trivial defects.

10.5. Purchaser shall provide Siemens with both on-site and remote access to the Products to the extent legally permitted and technically feasible. The remote access shall be provided through Purchaser's Internet Connection as is reasonably necessary for Siemens to provide services under this Contract, such as rectification of defects. The Purchaser agrees to allow connection to Siemens' remote service system. In the case of remote services, Siemens shall be entitled to access the data available with the Purchaser for the Products serviced by remote service.

10.6. Deliveries that exhibit material defects or lack an explicitly warranted characteristic at the time of risk transfer shall, at our discretion and within the normal processing times, be either repaired free of charge or replaced with non-defective deliveries. We shall assume ownership of any defective parts replaced. In the event of defects in the software supplied to Purchaser in connection with the Products – substantial non-conformance with the relevant software specification- Siemens provides Purchaser with an updated version of the software product eliminating such defect, which is or becomes available to Siemens. In the event modified standard software or developed individual software is defective, Siemens will until the provision of an updated version of the software eliminating the defect provide Purchaser with an acceptable by-pass or error eliminating solution, if this is possible with reasonable expenditure and if otherwise the commercial use of the Product is substantially affected.

10.7. Unless otherwise expressly agreed, we shall deliver repaired or replacement parts at our discretion FOB German port or FCA German airport in accordance with Incoterms 2010. We shall determine the port or the German airport to be used. We accept no responsibility for sending back or storing deliveries returned to us without our prior consent. For quality reasons, goods may be returned or taken back only in agreement with our order processing department.

10.8. If repair or replacement delivery is unsuccessful or if we allow a reasonable extension period to expire without rectifying the defect, the Purchaser may demand a price reduction or, if the defective delivery or the part of the delivery that is free from defects is of no use to the Purchaser, may withdraw from the contract.

10.9. No new warranty period shall commence as a result of repair, replacement or substitute performance/delivery.

10.10. When defects are reported, payments by the Purchaser may be withheld to an extent reasonably proportionate to the defects identified. However, the Purchaser may withhold payment only if a justifiable complaint has been brought. If the complaint made is unfounded, we are entitled to demand reimbursement from the Purchaser of our costs incurred thereby.

10.11. Further claims and rights of the Purchaser beyond those provided in this Item 10 asserted against us due to material defects or the lack of warranted features are excluded.

16. Contacts for reclamation

In case of any questions or problems with handling the system, please use the following contacts:

Contacts for questions and reclamations in Russia:

Siemens Healthcare LLC
Dubinskaya street 96
115093 Moscow,
Russia
Tel.: +7 495 737 12 52
Fax: +7 495 737 13 20

Customer support center
8 800 550 50 22
(24 hour free line)
service.ru@siemens-healthineers.com

Contacts for questions and reclamations in the country of Legal manufacturer:

Siemens Healthcare GmbH
Henkestr. 127,
91052 Erlangen
Germany
Tel.: +49 9131 84-0
Web: www.siemens.com/healthcare

17. Disposal of the system

CAUTION

Improper disposal of the system or parts of the system!

Pollution of the environment.

- ✦ System components hazardous to persons or the environment must be disposed of with care and in compliance with legally binding ordinances.
- ✦ Examples of environmentally relevant components are accumulators and batteries, transformers, capacitors, monitor picture tubes, and phantoms.
- ✦ For details contact your local customer service representative or your Siemens regional office.

Deed No. N 1380 /2021 Fr

I, Marcel Neumair, notary in Forchheim, notarize that the document overleaf was signed in my presence by

Mr. Ralf Rüdiger Körner,
born on September 24th, 1957,
with business address in D-91301 Forchheim, Siemensstraße 1,
personally known to me, notary,

acting according to his declaration for
Siemens Healthcare GmbH.

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Monika Schulz HC DI CT QT PQM Copyright © Siemens Healthcare GmbH 2018, 2018. Все права защищены. Конфиденциально.
Copyright © Siemens Healthcare GmbH 2018, 2018. Все права защищены. Конфиденциально.
Copyright © Siemens Healthcare GmbH 2018. 2018. Все права защищены. Конфиденциально.