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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality piece of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the KARL STORZ CALCULASE III. All the necessary details and actions are clearly explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective binder and keep it available for ready reference in a convenient and conspicuous location near the equipment.

All product illustrations, product descriptions and texts are the intellectual property of KARL STORZ SE & Co. KG. Their use and reproduction by third parties require the express approval of KARL STORZ SE & Co. KG.

All rights reserved.

It is recommended that the suitability of the products for the intended procedure be checked prior to the planned intervention.

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




**2 Controls, displays,
connectors and their uses**

- ① Keyswitch (0, 1)
- ② Emergency stop switch
- ③ Touchscreen
- ④ Ready/Standby button
- ⑤ Connection socket for footswitch
- ⑥ Fiber connection port/laser outlet aperture
- ⑦ Laser Active LED
- ⑧ Power cord socket
- ⑨ Power switch for turning the unit on and off
- ⑩ Potential equalization connector
- ⑪ Service port
(‘concealed’ – for service purposes only)
- ⑫ Interlock
(‘countersunk’ – for installation by service engineer; see Section 5.14)
- ⑬ Cooling system ventilation
- ⑭ Cooling water inlet
- ⑮ Cooling water fill level indicator
- ⑯ Cooling water outlet
- ⑰ Fiber tip
- ⑱ Protective cover around the laser fiber
- ⑲ Fiber connector
- ⑳ Footswitch with guard

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3 Symbols employed



3.1 Symbols on the unit

	Follow instructions for use
0/1	OFF/ON
	Laser emergency stop
	Ready/Standby button
	Fiber laser beam (laser fiber)
	Degree of protection against electric shocks: Applied part type B
	Footswitch connection
	Potential equalization connector
	Interlock connection
	Service port
	Vent
	Fill
	Empty
	This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	Manufacturer
	Electronic information product pollution control (China RoHS)

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3.2 User interface symbols and display elements

	Accessing procedure ①
	Procedure name ②
	Accessing menu (settings) ③
	'Info line' (error/warning/info) ④
	Pulse mode – 'short' ⑤
	Pulse mode – 'long' ⑤
	Pulse mode – 'burst' ⑤
	Indicator overview (totals): Energy – Pulse – Active time ⑥
	Reset 'totals' ⑦
	Energy indicator ⑧
	Frequency indicator ⑨
	Reduce/increase value
	Power indicator ⑩
	Fiber diameter indicator Number of applications ⑪
	Standby <=> Ready <=> Active ⑫
	NOTE: In the 'Active' status, the bar will be shown completely in yellow.
	Accessing previous setting level
	Accessing work screen

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3.3 Label and packaging symbols

	Manufacturer
	Catalogue number
	Temperature limit -10 °C ... +60 °C (storage temperature)
	Keep dry
	Fragile, handle with care
	Humidity limitation Storage 5%...95%
	Atmospheric pressure limitation Storage 500 hPa...1080 hPa
	Do not use if package is damaged
	Keep away from sunlight
	Consult instructions for use
	Federal (USA) law restricts this device to sale by or on the order of a physician.
	CE mark. The CE mark is a marking required for certain marketable industrial products in accordance with EU law. It is composed of the CE logo and, where applicable, the four-figure identification number of the responsible testing center (when the product belongs to a certain risk category).

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1 Images of the unit	V
2 Controls, displays, connectors and their uses	VII
3 Symbols employed	VIII
3.1 Symbols on the unit	VIII
3.2 User interface symbols and display elements	IX
3.3 Label and packaging symbols	X
4 System description	3
4.1 General information	3
5 Safety instructions	5
5.1 Explanation of Warnings and Cautions	5
5.2 Safety precautions at the site of installation	10
5.3 Safety precautions when operating the system	10
5.4 Intended purpose	12
5.5 Indications	12
5.6 Contraindications	13
5.7 Application risks and undesirable side effects	13
5.8 Intended use	14
5.8.1 User qualification	15
5.8.2 User profile	15
5.8.3 Patient profile	15
5.8.4 Operation	15
5.9 Operator's responsibility	16
5.9.1 Organizational measures	17
5.9.2 Before operating for the first time	17
5.9.3 Laser area	18
5.9.4 Warning signs for the laser area	18
5.10 Personal safety	19
5.11 Patient safety	21
5.12 Safety features of the system	21
5.13 Keyswitch	21
5.14 Emergency stop switch	22
5.15 Interlock	22
5.16 Fault detection	23
6 Installation and operating instructions	24
6.1 Unpacking and inspection	24
6.2 Basic equipment	24
6.3 Installing and connecting the unit	25
6.3.1 Filling the cooling system with cooling water and draining	25
6.3.2 Connecting the potential equalization	27

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6.3.3	Connecting the power cable	27
6.3.4	Switching on the main switch	28
6.3.5	Connecting the footswitch	28
6.4	Initial operation	29
6.4.1	Emergency stop switch	29
6.4.2	Switching the unit on	29
6.4.3	Starting the unit for the first time	30
6.4.4	Function test (with thermal paper)	31
6.4.5	Procedure list and procedure management	33
6.4.6	Menu	34
6.4.7	Changing to ready mode	40
6.4.8	Activating the laser beam	40
6.4.9	Removing the fiber	45
6.4.10	Switching the unit off	45
6.4.11	Trimming the fiber	46
7	Maintenance	47
7.1	Replacing the protective glass	47
7.2	Reprocessing	48
7.2.1	General warnings	48
7.2.2	Manual wipe-down disinfection	49
7.3	Maintenance and safety inspection	50
7.3.1	Maintenance	50
7.3.2	Safety inspection	50
7.4	Disposal	52
7.5	Repairing devices	53
7.6	Infection prevention	53
7.7	Limitation of liability	54
7.8	Warranty	54
8	Technical description	55
8.1	General	55
8.2	Troubleshooting	55
8.3	Technical data	59
8.3.1	Standard compliance	60
8.3.2	Directive compliance	60
8.4	Technical documentation	61
9	Replacement parts/accessories	63
10	Electromagnetic Compatibility (EMC) information	65
11	Subsidiaries	82

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4 System description

4.1 General information

By applying laser energy, lasers are used in the destruction of calculi and in tissue ablation, tissue vaporization and coagulation. The energy generated in the units is transferred through a probe/fiber to the target object, where it leads to coagulation, ablation/vaporization or fragmentation due to the absorption of laser energy. The fiber, in turn, is introduced into the surgical field through endoscopic instruments.

The CALCULASE III is a laser system that can be used for the cutting, coagulation, vaporization and ablation of soft biological tissue and the lithotripsy of calculi.

The device is operated and monitored via a touch display.

The laser has three status modes:

- Standby (settings can be adjusted)
- Ready (display is ready for operation; the footswitch is connected)
- Active (laser power is emitted as soon as the footswitch is depressed)

The activation of laser energy (Active mode) is performed via a footswitch, and is only possible in Ready mode. Standby mode is used to adjust the laser parameters and to display the various options.

Setting options for the CALCULASE III

- The laser has a maximum output energy of 35 Watt.
- Energy in the range of 0.2 to 4 Joule
- Frequency in the range of 4 to 30 Hz
- Three different pulse modes:
 - o Long pulse duration (450 µs ... 2000 µs)
 - o Short pulse duration (90 µs ... 450 µs)
 - o Burst mode (multiple pulses emitted in quick succession as a block, repeated with the specified frequency)

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There are eight preset procedures:

- Lithotripsy of a hard calculus (1.4 J; 6 Hz; short pulse duration) – standard setting following a system restart
 - Lithotripsy of a soft calculus (0.8 J; 6 Hz; short pulse duration)
 - Lithotripsy with a short recoil (1 J; 6 Hz; burst pulse mode)
 - Lithotripsy dusting (0.6 J; 15 Hz; long pulse duration)
 - Ablation (2 J; 10 Hz; short pulse duration)
 - Coagulation (1 J; 15 Hz; long pulse duration)
- Additional storage space is available for storing individual procedures.
- Sialendoscopy: 0.6 J; 4 Hz, burst pulse mode
 - Popcorning: 0.8 J; 15 Hz; long pulse duration

ⓘ NOTE: Simultaneous activation of the CALCULASE III with an HF surgical device is not envisaged.

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5 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

5.1 Explanation of Warnings and Cautions

Please read this manual carefully and follow its instructions precisely. The words Warning, Caution and Note convey special meanings. Wherever they are used in the instruction manual, please read the text that follows carefully, in order to ensure safe and efficient operation of the device. To make the signal words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician is at risk.

Disregarding a Warning or improper use of the laser could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

NOTE: A Note indicates special information about operating the device or clarifies important information.

WARNING: Risk of injury and damage to the products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals carefully and always follow the instructions given precisely. Check the functioning of the products used in combination.

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WARNING: The instruction manual must be kept with the CALCULASE III at all times. In addition to the instruction manual, the generally applicable and local rules on accident prevention and environmental protection must be available and must be observed.

All the safety and hazard information on the unit must be kept in a legible condition and replaced if necessary.

WARNING: The KARL STORZ CALCULASE III is a Class 4 laser according to IEC 60825-1. The laser may only be operated in designated laser areas on the authority of a surgeon.

WARNING: The laser may only be operated after receiving instruction by KARL STORZ or an authorized representative on the basis of this instruction manual. Further instructions to be followed before operating the unit for the first time are given on p. 17.

WARNING: The instruction manuals and interface specifications for medical devices and/or system components used in combination must be observed precisely.

WARNING: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective manuals or
- the intended use and interface specification of the devices used in combination permit this (cf. Clause 16 of IEC 60601-1-1).

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Safety instructions Warnings and Cautions

-  **WARNING:** Improper use which does not comply with the contents of the instruction manual may result in dangerous radiation.
-  **WARNING:** Do not open the device! Danger of electric shock. Have servicing carried out by authorized personnel only.
-  **WARNING:** Grounding reliability can only be achieved when the equipment is connected to a 'Hospital Only' or 'Hospital Grade' outlet (i.e., approved for use in an operating room environment). Routinely inspect the electrical plug and cord. Do not use if the inspection reveals damage.
-  **WARNING:** Inadvertent activation of the footswitch in READY mode can trigger laser pulses. For safety reasons, the laser device should be switched from READY mode to STANDBY mode between treatments.
-  **WARNING:** The electrical installations in the operating room in which the unit is connected and operated must comply with the applicable IEC standards or equivalent national regulations.
-  **WARNING:** Keep the unit out of reach of patients and do not touch the device and patient at the same time.
-  **WARNING:** Always unplug the unit before performing any maintenance work on it.
-  **WARNING:** Defective items of equipment must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.
-  **WARNING:** Caution – Use of the operator controls or setting options in ways other than those described here can result in dangerous irradiation.

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Safety instructions
Warnings and Cautions

-  **CAUTION:** Avoid allowing liquids to enter the unit. Do not store liquids on or above the unit.
-  **CAUTION:** Only operate the CALCULASE III with the voltage stated on the rating plate.
-  **CAUTION:** A laser unit, just as any other high performance medical device, requires special knowledge and care with regards to handling and use. The laser may only be operated by physicians who have been thoroughly trained in handling the device in accordance with this operating manual and are familiar with its therapeutic effect and possible risks.

 **NOTE:** Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

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The KARL STORZ CALCULASE III has been specially designed to minimize accidental exposure to hazardous laser radiation.

The laser unit is also equipped with a visible green aiming beam with a wavelength of 532 nm, which helps the surgeon to aim the laser beam.

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Reflected and direct eye exposure hazard

During the procedure, the laser emits high radiation power in the invisible spectral range of 2080 nm, which is hazardous to the human eye and skin.

As a precaution against accidental exposure to the output laser beam or its reflections, all personnel and patients must use appropriate safety eyewear (see Section 5.9 'Personal safety').

Do not stare into the laser beam or allow it to be reflected from any surface. Even rough metal can reflect the laser beam.

! WARNING: Never stare into the tip of the fiber. To prevent skin burns, do not aim the laser beam directly at unprotected skin.

! WARNING: Never stare directly into the aiming laser beam. Also ensure that the aiming beam is not directed at anyone's eyes. Although the beam is low powered, direct exposure can be hazardous to the eyes.

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
5.2 Safety precautions at the site of installation

The system may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

It is not intended for use in potentially explosive atmospheres. This means, among other things:

When using easily combustible and explosive inhalation anesthetics or mixtures thereof, the system must not be operated inside the demarcated hazard zone (Fig. 1). This also applies for easily combustible and explosive chemicals, e.g., skin disinfectants and fast-acting surface disinfectants.

The unit is equipped with a connector for potential equalization. It should be connected in accordance with the applicable national regulations.

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5.3 Safety precautions when operating the system

It is the user's responsibility to make sure the system is safe and operates properly before using it.

During treatment with the CALCULASE III system, the patient must be treated and kept under observation with the usual medical care. This includes constant monitoring of the progress of treatment, the patient's vital signs and anesthesia.

Any treatment must only be performed if visual observation of the action of the unit is ensured.

This unit has been tested and complies with the EMC limit values in accordance with the Medical Device Directive 93/42/EEC (CISPR 11 Class B and EN 60601-1-2). Please read the Electromagnetic Compatibility information in Section 10 on p. 65-78).

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The laser radiation of the CALCULASE III has a max. power of 35 W and a wavelength of 2080 nm, which is invisible to the human eye (the green aiming beam operates on a wavelength of 532 nm and is visible to the human eye). Incorrect use may lead to third degree burns even if the irradiation is unfocused.

Flammable drapes, surgical gowns, gauze and other ignitable materials must be kept clear of the laser beam path.

The use of non-flammable materials and instruments is advised. Flame retardant surgical drapes, gowns, etc., are recommended.

WARNING: There is a risk of fire or explosion if the LASER BEAM is used in the presence of flammable materials, solutions or gases or in an oxygenated environment. Some materials which are saturated with oxygen, for example, cotton, may ignite when exposed to the high temperatures produced during the INTENDED USE of the laser device. The solvents in adhesives and flammable solutions, which are used for cleaning and disinfection, should be allowed sufficient time to evaporate before the laser equipment is used. Please note that endogenous gases are flammable.

WARNING: Laser smoke may contain viable tissue particles. Smoke evacuation is necessary for applications that produce smoke gas.

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5.4 Intended purpose**Laser units**

Laser units are intended to provide laser radiation for cutting, coagulation, vaporization and ablation of biological tissue, as well as for lithotripsy of stones during surgical procedures. Laser units do not have body contact.

Laser fibers

Laser fibers, in combination with Calculase, Calculase II or Calculase III units, are intended for the transmission of laser radiation during surgical procedures. Laser fibers are surgically invasive or invasive (natural body orifice) and meant for transient use.

Accessories

Accessories are intended to support the laser units, or laser fibers. Accessories do not have body contact.

5.5 Indications

The use of the laser units and laser fibers with the corresponding accessories is indicated when laser surgery is indicated from the medical view of the treating physician, in the treatment of the following conditions:

- benign prostatic hyperplasia (BPH) / benign prostatic enlargement
- bladder tumors
- urolithiasis
- cholelithiasis
- sialolithiasis
- condyloma
- stenosis in urinary tract
- bladder outlet obstruction (BOO)
- upper tract urothelial tumors
- removal of eroded sling / mesh in urinary tract

The laser units with the corresponding accessories are suitable for application parts, which require the energy in the 'Technical Data' chapter of the instructions for use for successful performance.

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5.6 Contraindications

The use of the laser units and laser fibers with the corresponding accessories is contraindicated if, from the medical view of the treating physician, the surgical method is contraindicated, or if the patient's general condition does not allow for surgery or anesthesia.

The use is contraindicated if, from the medical view of the treating physician, the laser unit and laser fibers with the corresponding accessories are not compatible with the successful performance of the planned intervention due to their technical design.

The laser units are not used in body contact but provide energy and power for the corresponding application parts, which are transferred to the tissue/stones via the laser fibers.

The laser units and laser fibers with the corresponding accessories must not be used for interventions in direct contact with the CNS (central nervous system) and central circulatory system. Furthermore, intravascular use is contraindicated.

5.7 Application risks and undesirable side effects

The following side effects / complications have been recorded in the literature in relation to interventions where a Holmium:YAG laser was used:

- acute urine retention
- bleeding requiring coagulation
- bladder contracture / bladder neck stricture
- bladder perforation
- bladder irritation
- bladder tamponade / clot retention
- blood transfusion
- dysuria
- auxiliary procedures
- bile leakage
- haemobilia
- urinary tract infection
- incontinence/ stress urinary incontinence

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- capsular perforation
- cardiac events
- meatus stenosis
- perioperative mortality
- postoperative gross hematuria
- postoperative fever
- re-hospitalization
- re-catheterization
- reoperation
- storage symptoms
- retrograde ejaculation
- recurrence
- stone migration
- strictures in the urinary tract (urethra, ureter)
- ureter perforation
- urethral perforation
- bladder mucosal injury (superficial)
- ureteral orifice injury
- reduced erectile function

In addition, clinical literature has shown an increased rate of side effects for patients with abnormalities of blood coagulation.

5.8 Intended use

The unit must only be operated with accessories which have been designated by KARL STORZ as suitable for the unit or the safe technical functioning of which has been proven.

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5.8.1 User qualification

The CALCULASE III must only be used by physicians and medical assistants with the corresponding professional qualifications who have been instructed in using the device. Instruction is given by authorized persons. The local KARL STORZ representative is available to provide you with more information on instruction.

5.8.2 User profile

- Recognized medical skills in the field of application
- Sufficient understanding for the rational assessment of the medical intervention and use of the device.
- Familiarity with technical devices
- Suitable command of the language used by the device and in the instruction manual
- Knowledge of the contents of the instruction manual
- Adequate knowledge of the purpose and whereabouts of the medical products
- No physical impairment which makes it difficult or impossible to perceive acoustic or visual alarm signals

5.8.3 Patient profile

The patient population is not limited. The medical device can be used for all patients requiring laser procedures.

5.8.4 Operation

The unit must only be used in hospitals and by physicians in the ambient conditions specified in the technical data.

Other conditions of use:

- Frequency of use: one or more times per day
- Duration of use: from a few minutes to several hours per day
- Positioning: on a level, vibration-free surface
- Can only be moved by being installed on trolley or ceiling supply units

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- Combination: can be used on the patient at the same time as other devices required for the operation
- Lasers are not designed for use in potentially explosive atmospheres. That means that the system cannot be operated in fields of application where flammable or explosive inhalation anesthetics, chemicals (e.g., skin disinfectants and rapid-acting surface disinfectants) or mixtures thereof are used.
- In order to prevent damage to the fiber system, laser fibers are not intended for use with cutting systems such as razors or morcellators.

5.9 Operator's responsibility

The operator and safety officer have responsibility for carrying out all the safety measures to prevent any risk to the patient, attending physician and other people present during use of the laser.

In addition, the national accident prevention rules and regulations applicable to the place of use must be adhered to.

In Germany, the valid accident prevention regulations (German Social Accident Insurance Regulation 11), the OStrV [German regulation on artificial optical radiation] and the associated Technical Rules under the German OSH Ordinance on Artificial Optical Radiation (TROS laser radiation) must be observed.

The operator must keep the following documents accessible at all times and must update them regularly (Germany only):

- Inventor list
- Equipment book and/or equipment chart
- Instruction manual

Accidents and personal injury must be reported to the responsible authorities without delay.

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5.9.1 Organizational measures

The required personal protection equipment must be provided by the operator.

All the safety devices fitted must be checked at regular intervals.

i NOTE: Secure the unit against unauthorized use by pulling the key out of the keyswitch. The key can only be removed when in the '0' position.

5.9.2 Before operating for the first time

The initial installation of the CALCULASE III must be carried out by a trained, authorized person. The system must first be filled with cooling water, otherwise start-up will not be possible (see also the caution note on page 24).

The guidelines of the country of use for installing and setting up laser systems must always be followed.

The following points apply to Germany:

- The unit must be inspected and approved by the procurement office.
- The procurement office or operator must appoint a laser protection officer in writing.
- The Employer's Liability Insurance Association and the health and safety authority must be notified of the laser unit.
- All the structural safety requirements in the laser area must be met.
- The unit must be handed over ready for operation and the staff trained by KARL STORZ or an authorized representative.

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5.9.3 Laser area

The laser area is the area in which the levels for maximum permissible irradiation may be exceeded. This means the entire operating room unless building measures have been carried out to limit the area. The possibility of unintentional deflection of the laser beam must be taken into account. The access routes to the laser areas must be identified with warning lamps. A flashing light is recommended.

Signs must be placed in the laser area giving an explicit warning about the intensive light and the laser beam used.

5.9.4 Warning signs for the laser area

Position the warning signs conspicuously.

Only allow access to the CALCULASE III laser area to those members of staff who are absolutely necessary for the treatment and who have been thoroughly trained in the necessary safety measures.

CAUTION: In the laser area, no reflective objects, e.g., surgical instruments, mirrors, watches or jewelry are permitted which might reflect the laser.

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5.10 Personal safety

- Ensure that all members of staff in the laser area are familiar with controlling the CALCULASE III and know how to shut down the system rapidly in an emergency.
- Never look directly into the light beam coming from the laser fiber, even if you are wearing safety goggles.
- Never direct the laser fiber around the room.
- When the CALCULASE III is not being used, switch it off immediately and remove the key.

Laser safety goggles

Laser safety goggles prevent the eyes from being harmed by intensive light during a light pulse. The goggles are equipped with side guards for additional safety against lateral reflections.

WARNING: During treatment, the physician and everyone who is present in the laser area must wear special laser safety goggles for a 2080 nm wavelength (e.g., laser safety goggles 27 750095). The laser safety goggles for the CALCULASE III laser system must meet at least the following conditions:
I 2080 nm LB3 (in accordance with EN 207),
Marking declaration:
I: Laser mode (pulse laser)
2080 nm: wavelength at which the laser safety goggles provide protection.
LB3: protection level

CAUTION: Non-compliance can lead to irreversible eye damage. Please verify the perfect condition of the laser safety goggles before use – the goggles must be free of any mechanical damage. Only the type of protective eyewear for the indicated wave length spectrum may be used.

NOTE: The frame or arms of the goggles should have a label indicating the specific wave length spectrum they are designed for. The wavelength on the label must correspond to or include the wavelength of the laser.

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5.11 Patient safety

The safety of the patients is ensured primarily by well-trained staff and a laser area equipped according to the regulations. The patients must be informed about the safety measures which are important for them and about the type of treatment.

- Ensure that the patient's eyes are protected against the laser beam.
- If the patient is not able to wear safety goggles, provide him or her with opaque eye protection that entirely prevents the laser beam from reaching his or her eyes.

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5.12 Safety features of the system

The safety standards IEC 60601 and IEC 60825-1 determine the general requirements relating to the design and manufacture of the KARL STORZ CALCULASE III.

In compliance with international standards, the CALCULASE III is equipped with the following:

- Energy indicator
- Keyswitch ①
- Emergency stop switch ②
- Interlock connection ③

Proper signs and labels (according to IEC 60825), see Figures on the left (Fig. 3)

! WARNING: The unit is grounded by a grounded conductor in the power cable. This grounding is necessary for safe operation.

Before switching on the unit, all the safety devices must be properly fitted and fully functional.

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
5.13 Keyswitch

The keyswitch ① is used to turn the unit on and off (Fig. 4). The system can only be operated with the key provided by KARL STORZ. This prevents unauthorized use of the unit. Always remove the key after switching the unit off and store it in a safe place.

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5.14 Emergency stop switch


The red round button is used for shutting down the unit rapidly in an emergency. If you press the emergency stop switch ②, the laser emission will be disconnected immediately (Fig. 5).

Turn the unit to pull out the emergency stop switch ② again before switching the unit back on. When the emergency stop switch is pressed, the system can be switched on but laser activation is not possible.

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5.15 Interlock

For extra safety, the CALCULASE III is equipped with a connection point for a laser radiation interlock. This connection point ⑩ is located on the back of the unit (Fig. 6). An external switch is connected to this outlet, which enables remote interlocking.

The external switch is connected to the door of the operating room for example, so that the laser beam is interrupted when the door is opened. The unit then returns to standby mode and the error message 'interlock' is output.

To resume operation, the external switch contact must be closed and the unit turned to READY mode once again. Press the footswitch to reactivate laser radiation.

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

NOTE: The interlock must only be installed by authorized specialist personnel.

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5.16 Fault detection

The unit is equipped with a self-check routine that continuously monitors system operation. If a 'fault condition' is detected:

- an error message or warning message (info message) is shown on the screen (see figures).
- a signal tone is emitted.
- the system ceases laser emission to ensure user and patient safety and to prevent damage to the unit.

A list of error messages, their meanings and possible remedies can be found in Section 8.2.

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6 Installation and operating instructions

6.1 Unpacking and inspection

Carefully remove the device and the accessories from the packaging. Check for missing items and for evidence of shipping damage.

Please file any complaints with the manufacturer or supplier immediately.

If possible, retain the original packing materials for later use; these can be used if the device has to be transported.

CAUTION: The initial installation of the CALCULASE III must be carried out by qualified personnel. To ensure safe operation, the CALCULASE III should be filled with cooling water.

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6.2 Basic equipment

The CALCULASE III basic equipment comprises:

UL 300	CALCULASE III
ET21-1002	Power cable
27 7501 24	Footswitch
27 7502 35	Set of keys for the CALCULASE III
UL 101	Interlock bridge (the interlock is delivered inserted into the interlock connection ⑫ on the back of the unit as standard)
27 7500 95	Laser safety goggles Ho:YAG laser 2080 nm
UL 102	Filling set for cooling water
96136029 DER	Instruction manual

NOTE: Optional accessories can be found under 'Recommended Accessories'.

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

6.3 Installing and connecting the unit

Set up the CALCULASE III on a flat surface as a free-standing unit (ideal for a KARL STORZ equipment trolley) (Fig. 10).

NOTE: The CALCULASE III may only be used in medical environments set up in accordance with applicable national regulations for laser use.

NOTE: The KARL STORZ CALCULASE III is an extremely mobile system thanks to the desk-top design. The laser can be transported by its special foldaway handle (Fig. 11).

6.3.1 Filling the cooling system with cooling water and draining Filling the cooling system

When filling the cooling system, please proceed as follows:

1. Attach the supplied hose without funnel to the top left-hand connection (13) (Figs. 12, 13) (venting).
2. Attach the supplied hose with funnel to the top right-hand connection (14) (Figs. 12, 13).
3. Pour distilled or deionized water (approx. 1.5 liters) into the funnel.

NOTE: The laser crystal that is flushed with the cooling water may be damaged as a result of the ions in the water. For this reason, the device is already factory-fitted with a deionization filter (ion exchanger) that must be replaced once a year as part of the safety inspection.

4. The fill level is shown at display (15) next to the connections (Fig. 12). Fill the cooling system up to the top marking.

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5. Disconnect the hoses and store in a dust-free environment for the next safety inspection.
6. Start the unit to check whether the fill level is sufficient. If necessary, switch off the unit, top up with more cooling water and start again.

Draining the cooling system

Drain the water from the cooling system for maintenance or transport.

- i NOTE:** During transport or storage, there is a risk that the cooling water may freeze and cause damage to the unit. Therefore, the cooling water must be drained prior to transport or prolonged periods of storage.

To do so, proceed as follows (Fig. 14):

1. Attach the vent hose to the top left-hand connection 13.
2. Connect one of the two bottom right-hand drainage connections 16 to the hose piece without funnel.
3. Guide the end of the hose to a sufficiently large container.
4. Allow the cooling water to completely drain out of the CALCULASE III. Tip the unit onto its two rear feet to remove any residual water.
5. Connect the hose to the other bottom right-hand connections and drain as described under points 2 and 3.
6. Remove the vent hose and the drainage hose.

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
6.3.2 Connecting the potential equalization

The CALCULASE III is equipped with a connector for potential equalization ⑩. The potential equalization is the conductive electrical connection of device housings both between themselves and to ground. It ensures that the devices always contain the same potential, even in the event of electrical faults. Potential equalization is prescribed for certain operating rooms.

The device's ground line should be installed by a qualified electrician (Fig. 15).

Protection Class I device

WARNING: To avoid the risk of an electric shock, this device may only be connected to a power supply network with a protective conductor.

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
6.3.3 Connecting the power cable

Insert the power cable into the power socket ⑧ to the stop (Fig. 16).

Only operate the CALCULASE III with the voltage stated on the rating plate (Fig. 17).

The power cord connector must be pulled out to ensure an all-pole disconnection from the mains power supply. For this reason, the unit must not be set up in a manner which makes it difficult to disconnect it from the mains power supply.

NOTE: If the laser is operated with low voltage (< 120 V), it will not be able to deliver the full output energy (see parameter matrix).

WARNING: The unit must not be operated in potentially explosive atmospheres.

CAUTION: Do not connect the laser either to a multiple power socket outlet or a socket that can be switched off.

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


6.3.4 Switching on the main switch

After connecting the power cable to the power supply, switch on the main switch ⑨ (Fig. 18).

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6.3.5 Connecting the footswitch

Insert the footswitch connecting cable into socket ⑤ (in alignment with red dot) (Fig. 19).

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6.4 Initial operation

6.4.1 Emergency stop switch

To begin laser operation, the emergency stop switch (2) must be unlocked by turning it (Fig. 20).

NOTE: When the emergency stop switch is pressed, a corresponding note will appear in the GUI note cell.

Press the emergency stop switch (2) in an emergency; the laser is switched off instantly. To resume operation, the switch must be unlocked.

WARNING: The keyswitch (1) and the power switch (9) must be used for normal system shutdown. Pressing the emergency stop switch (2) immediately stops laser emission. If the emergency stop switch (2) has been pressed, it will need to be unlocked again by turning.

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6.4.3 Starting the unit for the first time

If the CALCULASE III is being started for the first time, the screens for setting the language and the time and date will appear in succession.

See the 'Menu' section starting on page 33 for information on setting the language and time/date.

After the unit has been switched on, a white screen will appear for a few seconds, followed by the screen with the KARL STORZ logo as shown on the left (see Fig. 22).

NOTE: A signal tone will be emitted shortly before the end of the switch-on process. (check whether a signal tone is emitted).

CAUTION: In order not to damage the sensitive touchscreen surface, do not use any objects to tap the screen.

NOTE: Following initial installation, the unit will start up each time with the standard Lithotripsy of a 'hard calculus' setting, with the following parameters: 6 Hz, 0.8 J and short pulse duration. Default settings cannot be overwritten.

Menu structure

The menu structure of the CALCULASE III is split into 3 areas:

- Working area
- Procedure list and procedure management (see Section 6.4.5)
- Settings (see Section 6.4.6 starting on page 33)

The working area can have the status of Standby (see Figure at the bottom left), Ready or Active. If the unit's protection systems detect an error, the system will return to Standby status.

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6.4.4 Function test (with thermal paper)

1. Make sure that the laser has been correctly filled with water and that the footswitch is connected (see Sections 6.3.1 and 6.3.5).
2. Make sure that the emergency stop switch is unlocked and that the key is in the keyswitch (see Sections 6.4.1 and 6.4.2).
3. Turn the key clockwise to position 1 (see Section 6.4.2).

! WARNING: The fiber produces an aiming beam as soon as it is switched on. Never look directly into the laser beam.

! CAUTION: The tip of the fiber must be free of dust and grease.

! CAUTION: The use of fibers with a defective or soiled connector at the coupling results in the reduced emission of laser energy and destruction of the protective glass which must then be replaced.

5. Aim the fiber at a sheet of thermal paper (Fig. 26).
6. At a setting of 6 Hz/0.8 J and with a short pulse mode, switch the device from Standby to Ready.
 - The tip of the fiber and the thermal paper should be approximately 10 cm apart; the aiming beam should be approximately 20 mm in diameter.
 - The aiming beam emitted should be round and homogenous.

! WARNING: Only aim at the front side of the thermal paper!

! WARNING: Thermal paper is flammable! Stop activation immediately once the blackening is the size of the aiming laser beam! Stop activation immediately if smoke develops!

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7. Press the footswitch to activate the laser.
 - There must be no visible effect on the thermal paper (Fig. 27).
8. Release the footswitch again.
9. Reduce the distance between the tip of the fiber and the thermal paper to approx. 5 cm; the pilot laser beam should now be approx. 8 mm in diameter.
 - The aiming beam emitted should be round and homogenous.
10. Press the footswitch to activate the laser.
 - The thermal paper is blackened after a few pulses (Fig. 28).
11. Change to a setting of 6 Hz/2.0 J.
 - The tip of the fiber and the thermal paper should be approximately 10 cm apart; the aiming beam should be approximately 20 mm in diameter.
 - The aiming beam emitted should be round and homogenous.
12. Press the footswitch to activate the laser.
 - The thermal paper is blackened after a few pulses (Fig. 29).
 - On successful completion of the function test, there must be a small spot of approx. 8 mm in diameter and a large spot of approx. 20 mm in diameter on the thermal paper (Fig. 30).

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

6.4.5 Procedure list and procedure management

Pressing the navigation button from the 'working area' opens a list of procedures for the corresponding screen. From here, new procedures can be created or existing ones deleted or edited.

There are six preset procedures (Fig. 31):

- Lithotripsy of a hard calculus (0.8 J; 6 Hz; short pulse duration) – standard setting following a system restart
- Lithotripsy of a soft calculus (0.4 J; 6 Hz; short pulse duration)
- Lithotripsy dusting (0.4 J; 25 Hz; short pulse duration)
- Lithotripsy with a short recoil (1 J; 6 Hz; burst pulse mode)
- Ablation (2 J; 10 Hz; short pulse duration)
- Coagulation (1 J; 15 Hz; long pulse duration)

These are factory settings and cannot be deleted.

The CALCULASE III has the option of saving values. The set values for power and frequency are saved within a so-called procedure, together with the underlying operating mode. A list of these procedures can be accessed via the navigation area in the work screens (Fig. 32).

Procedure list

The procedure list features the existing procedures that have already been created. Existing procedures in the list can be modified, or new procedures created.

The procedure list is always in chronological order, with the preset standard procedures always at the top of the list.

Procedure list editing mode

If a procedure is selected using the checkbox, the procedure list will switch to editing mode. The title bar will be replaced by an action bar (Fig. 33).

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6.4.6 Menu

To open the menu screen, tap on the  'Menu' symbol (see e.g., Section 3.2).

The Settings screen appears (see Fig. 35).

The Settings screen features the following sub-menus:

- Pilot laser
- Protocol
- Device and fiber information
- License information
- Menu > Settings
 - Language
 - Date and time
 - Brightness
 - Service
 - Sound volume

There are multiple levels that can be accessed within the menu area: the unprotected area (unit settings), the Settings area and the Service area (password-protected).

Pilot laser (Fig. 36)

The brightness of the aiming beam can be adjusted in this menu ('+' to increase and '-' to reduce) Setting range 'off', 1...10.

NOTE: If the aiming beam is switched off, it will be switched on again for safety reasons once the unit has been restarted.

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

Protocol (Fig. 37)

This screen contains an overview of the recorded application data (totals for energy, pulse, active time).

Device and fiber information

This screen contains information on the device and the fiber, e.g., software configuration, hardware version, operating hours and number of applications of the connected fiber.

License information

This screen contains information on the licenses.

Menu >

Additional sub-menus can be opened from here.

> Language

This sub-menu can be used to set the language to appear on the screen.

Tap the required language to select it (the currently selected language is shown in orange). 15 languages are currently available.

> Date and time

This menu is self-explanatory.

> Brightness (Fig. 38)

The desired brightness can be set here.

> Service (password required – Fig. 39)

After entering the final digit in the password, the system will automatically switch to Service mode. No separate password entry confirmation is required.

The Service menu is reserved for the authorized Service and, for this reason, access to this menu is password-protected.

> Sound volume (Fig. 40)

The volumes for warnings sounds, information sounds and keypad sounds are individually adjustable (within the range of 1-10) or can be muted with one exception (see following information).

NOTE: The laser activation sound cannot be permanently muted.

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The laser is equipped with a fiber recognition system (RFID technology) that provides a permissible parameter profile according to the fiber diameter being used. Information on the 'detected' fiber is shown at the bottom left of the screen.

Parameter sets are selected via the fiber connector coding (see 'Parameter matrix' tables starting on p. 39).

The integrated RFID technology allows for automatic identification and reprocessing checks (Fig. 41):

- The fibers are detected.
- The number of applications is shown.
- The minimum and maximum possible setting values are applied automatically.

After the unit is switched on (exception: initial installation), the system automatically adopts the standard Lithotripsy of a 'hard calculus' setting, with the following parameters: 0.8 J, 6 Hz and short pulse duration.

Connecting the fiber

The unit cannot be used until a fiber has been connected.

If no fiber is connected, the words 'no fiber connected' will be shown at the bottom left of the screen.

To connect the fiber, it will firstly need to be opened (Fig. 42), the lever on the connection socket ⑥ moved to the left or right (Fig. 43) and the fiber simultaneously turned to secure (Fig. 44). The unit can be operated as soon as the fiber has been connected.

WARNING: Only clean fibers in perfect condition that have been sterilized according to the instructions may be used for treatment. The reprocessing of the fibers must always be monitored.

CAUTION: In order to prevent damage to the unit, do not touch the proximal end of the fiber. If the fiber connector is soiled, this can lead to the destruction of the protective glass or even the destruction of the fiber coupling focusing lens.

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Fibers are identified by the color of the connector.

Fiber	Color
230 µm	blue
365 µm	green
600 µm	red

Only KARL STORZ CALCULASE III fibers (p. 63) may be used with this device. In particular, no electrically conductive fibers may be connected.

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The following tables provide an overview of the correlations between fiber diameter, ablation volume, repulsion, wear of the fiber tip, possible flushing and possible bending of a flexible uretero-renaloscope.

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Curve (when using a Flex X)			
without fiber	230 µm	365 µm	
270°	225°	180°	

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! WARNUNG: Do not activate the laser if the fiber is faulty.

! WARNUNG: Avoid extreme bending or tight winding of the fiber. This could cause the fiber to break. The minimum bending radius is between 1.5 and 3.5 cm, depending on the fiber (230 µm fiber: 1.5 cm, 365 µm fiber: 2.0 cm, 600 µm fiber: 3.5 cm). Ensure that this is also taken into consideration for storage.

! WARNUNG: Pushing the laser fibers forward in the instrument under poor visual conditions can result in the perforation of the bladder or ureter as a result of the mechanical force exerted.

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The fiber must be checked at the beginning of every laser treatment. To do so, the aiming beam must be rechecked and reset every time it is used. The brightness of the aiming beam can be adjusted in the 'Settings' menu under the 'Pilot laser' menu item ('+' to increase and '-' to reduce) Setting range 0 ('off') ... 10.

! WARNING: As the aiming beam follows the same path through the laser transmission system as the WORKING BEAM, it offers a good means of checking the integrity of the laser transmission system. If the AIMING DOT does not appear at the distal end of the laser transmission system, its intensity is too weak, or it appears diffuse, this is a possible sign that the laser transmission system is damaged or not functioning clearly.

1. Inspect the fiber sheath surfaces for damage and porous areas before each use.
2. Switch the laser equipment on.
3. Connect the laser fiber.
=> The aiming beam is ready for use.
4. Check the distal reproduction of the aiming beam light.
=> Optical laser fiber transmission is OK if the aiming beam produces a round, evenly bright mark at a distance of 5 to 10 cm from the distal fiber end (Fig. 45).
=> Optical laser fiber transmission is not OK if the mark produced is not round, is unevenly bright, or is extremely irregular (comet tail).
5. Check the distal laser fiber surface and trim the laser fiber if necessary.
6. Replace the laser fiber if treatment results prove inadequate.

The laser light must only emerge at the distal end of the fiber. If you can identify multiple light points, this means the fiber is damaged.

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WARNING: Never stare into the tip of the fiber. To prevent skin burns, do not aim the laser beam directly at unprotected skin.


CAUTION: Always check the aiming beam in the standby mode and not in the READY mode. The beam emitted by the laser fiber should be round and homogeneous.

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The frequency is set in the working area using the '+' and '-' keys in the frequency field. There are up to eleven different frequency levels to choose from depending on the fiber being used (Fig. 46).

NOTE: At low voltage (< 120 V), not all parameters can be adjusted.

NOTE: The frequency and energy settings are interdependent. When the energy (or frequency) is changed, the other value may change accordingly.

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
The energy is set using the '+' and '-' keys in the energy field (Fig. 47). Up to eleven different energy levels are available here depending on the fiber and specified frequency.


NOTE: At low voltage (< 120 V), not all parameters can be adjusted.

The laser is always in standby mode during and after the parameter adjustment. The word 'Standby' is at the bottom right of the work screen.

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6.4.8 Activating the laser beam

The unit is in the READY mode. Press the footswitch to activate the laser beam.

The activation/output of laser energy is displayed by the word 'Active' shown in yellow (Fig. 50), by the 'Laser Active LED' ⑦ (see page VI) and by an activation tone. The set energy and frequency are highlighted in yellow.

WARNING: For safety reasons, the laser should be switched to standby mode during extended periods of (foreseen) non-activation.

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

The CALCULASE III is equipped with a fiber recognition system. A permissible frequency/energy combination is provided depending on the fiber being used, from which the user can select the corresponding capacity (see also p. 36).

! **CAUTION:** The use of fibers not approved by KARL STORZ for use with the unit may damage the fiber recognition system and cause the unit to fail.

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Parameter matrix for fibers »230 µm«
(short pulse) – energy [J] against frequency [Hz]

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Parameter matrix for fibers '230 µm' (long
pulse) – energy [J] against frequency [Hz]

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Parameter matrix for fibers '230 µm' (burst
pulse) – energy [J] against frequency [Hz]

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Parameter matrix for fibers '365/600 µm'
(short pulse) – energy [J] against frequency [Hz]

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Parameter matrix for fibers '365/600 µm' (long
pulse) – energy [J] against frequency [Hz]

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
6.4.9 Removing the fiber

For safety reasons, switch the device to standby mode; to do so, press the ready/standby key ④ (see p. VI). Then remove the fiber.

Seal the reusable fiber with the cover cap for reprocessing (Fig. 51).

Inspect the reusable fiber for defects prior to reprocessing. Dispose of defective fibers or reprepare where possible (see CALCULASE III fibers instruction manual, Art. no. 97000238).

NOTE: You will find information on the sterilization of reusable fibers in the instruction leaflet accompanying the fibers.

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
6.4.10 Switching the unit off

If the unit is not yet in standby mode, press the Ready/Standby button ④ (Fig. 52) to switch the unit to standby mode.

The unit is switched on (position '1') and off (position '0') using the keyswitch ① on the front of the unit (Fig. 53).


NOTE: The unit is switched off when the keyswitch is in position '0'.

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i HINWEIS: Abisolierlänge 10-15 mm.

6.4.11 Trimming the fiber

To improve the quality of the beam, the fiber can be trimmed. First the fiber must be stripped.

The fiber is then inserted into the indicated wire strippers depending on its diameter, both arms pressed together, and the fiber pulled out of the closed wire strippers.

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
To break the fiber cleanly, it is scored with a ceramic scoring pen (27 7502 84). It should then be held taut over the silicon pad (27 7502 85) and broken.

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
7 Maintenance

7.1 Replacing the protective glass

Carry out this step with clean hands and the utmost care so as not to damage any optical parts of the unit.

The protective glass is replaced as follows:

- Switch off the unit and disconnect from the power supply.
- Open the side flap (Fig. 58).
- Remove the protective glass (Fig. 59) and check for visible damage (Fig. 60).
- If the glass is damaged, replace the protective glass including the metal carrier. When doing so, take care not to touch the protective glass.
- Close the flap.

① **NOTE:** If you require assistance, please contact your service representative.

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7.2 Reprocessing

7.2.1 General warnings

 **WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

 **WARNING:** When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

 **CAUTION:** When preparing and using solutions, follow the chemical manufacturer's instructions, paying close attention to proper concentration, exposure time and service life. Incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

 **WARNING:** National laws and regulations must be observed.

 **NOTE:** The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

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7.2.2 Manual wipe-down disinfection

WARNING: Before cleaning or disinfecting, unplug the power cable. Personnel are only allowed to clean the exterior surfaces. Never open the unit.

CAUTION: When cleaning the unit, make sure that no liquid penetrates into the unit. Take special care of the fiber connection socket/laser outlet aperture!

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications concerning material compatibility and exposure time must be observed. At the end of the necessary exposure time for the disinfectant, wipe the surface with a dry low-lint cloth. Do not use sprays because no liquid must be allowed to enter the interior of the unit.

Reprocessing of accessories

Extensive accompanying documentation is available for the accessories (Art. no. 97000238). It is essential that you refer to these documents for further information on handling and reprocessing.

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7.3 Maintenance and safety inspection

WARNING: At least once per year (within 12 months), maintenance and a safety inspection/repeat inspection of the complete laser system as defined by IEC 62353 must be carried out by a service technician authorized by KARL STORZ.

NOTE: Service and maintenance work must not be carried out when the equipment is being used on a patient.

7.3.1 Maintenance

Annual maintenance is essential. Regular maintenance can contribute to identifying potential problems before they become serious, as well as enhancing the device's reliability and extending its useful operating life.

Maintenance services can be requested from your local representative or from the manufacturer.

Maintenance of the CALCULASE III must include the following work at least:

- Changing the coolant
- Replacing the particle filter
- Replacing the ion exchanger

7.3.2 Safety inspection

WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, safety inspections/repeat inspections as defined by IEC 62353 must be carried out on this unit by a service technician authorized by KARL STORZ at least once a year.

As standard, information is shown in the display 60 days before a safety inspection is due (Fig. 61).

Please get in touch with your local KARL STORZ contact or Service to arrange an appointment for a safety inspection.

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One year after the last safety inspection, a message will be displayed indicating that a safety inspection is overdue (Fig. 62). Please get in touch with your local KARL STORZ contact or Service to arrange an appointment for a safety inspection, so that the system can continue to work optimally.

The inspection includes accessories, wearing parts, disposable articles and the interfaces connected to the unit that are used to operate the unit if this accessory equipment has any impact on the unit's operational safety.

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check safety-relevant labels and warning stickers for legibility.

Electric measurements

- Inspect the device safety fuses
- Measure protective ground resistance according to IEC 62353: For limit values, please refer to the current standard.
- Measure earth leakage current according to IEC 62353: For limit values, please refer to the current standard.
- Measure touch current according to IEC 62353: For limit values, please refer to the current standard.
- Measure patient leakage current according to IEC 62353: For limit values, please refer to the current standard.

Functional test

The test for proper operation must include the following tests at least:

- Emergency stop
- Interlock switch (remote interlock)
- Displays and control keys
- Checking the laser emission using thermal paper

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Documentation

The safety inspection/repeat inspection and the results must be recorded in the medical device logbook.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during the safety inspection/repeat inspection, the unit must not be operated until these defects or shortcomings have been eliminated by qualified technical servicing.

NOTE: For detailed information on the range and implementation of the safety inspections/ repeat inspections, please refer to the current version of the relevant Service Manual.

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7.4 Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the device as electronic scrap.

Please ask KARL STORZ SE & Co. KG, a KARL STORZ subsidiary, or your authorized dealer about the collection point in your area. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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7.5 Repairing devices

Only personnel authorized by the manufacturer may repair devices. Interventions described in these instructions for use are exempt from this rule.

In Germany, you can contact the following address for repairs:

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Take-Off Gewerbepark 83

78579 Neuhausen

Service hotline: +49 7461 708-980
E-mail: technicalsupport@karlstorz.com

In other countries, the respective KARL STORZ branches or specialist dealers are responsible.

7.6 Infection prevention

Contaminated devices may not be shipped. To prevent contact infections and airborne infections, devices must first be decontaminated. KARL STORZ reserves the right to send back contaminated devices.

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7.7 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the device is connected and operated meet the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

7.8 Warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see Chapter 'Subsidiaries'), even during the warranty period.

Opening of the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair or modification performed during the warranty period shall void all warranty.

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8 Technical description

8.1 General

In the event of the device malfunctioning, the troubleshooting list provided in this section will assist the user in finding the possible cause. The troubleshooting list indicates how to find and rectify the various malfunctions.

WARNING: High voltage is present in the device's interior. The device must only be opened by an authorized KARL STORZ Customer Service technician.

WARNING: Maintenance work that is not mentioned in this instruction manual must only be performed by technical personnel that are authorized by KARL STORZ.

WARNING: Maintenance work or modifications to the laser unit which are carried out contrary to this instruction may result in personal injury, damage to the unit and loss of any guarantee rights.

WARNING: Always unplug the unit before performing any maintenance work on it.

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8.2 Troubleshooting

The troubleshooting tables do not attempt to list all possible system failures. Any fault not listed, should be referred to KARL STORZ customer service.

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Possible notes and errors

Error code x: Contact Service.

Emergency stop switch depressed. To continue, turn the emergency stop switch counter-clockwise to release it.

Error: Connected fiber not recognized. Use a different fiber.

The fiber can no longer be used.

Last utilization cycle of the fiber. Discard after use.

Undetermined error x: If the error occurs again, contact Service.

Error: Energy setpoint not reached: Re-activate the laser. If the error persists: contact Service.

Error: Fiber connector overheating: Use a different fiber. If the error persists: Consult instruction manual.

Guaranteed light service life almost reached: Contact Service to arrange a light replacement.

Guaranteed light service life has been exceeded. Contact Service to prevent a laser failure.

Interlock active. Release interlock to continue.

Error in cooling system: Check water fill level and restart the unit.

No footswitch connected. Connect footswitch to continue.

Laser above operating temperature. Please wait until the system has cooled down to operating temperature.

Footswitch depressed in standby mode. Release the footswitch.

Error: Aiming beam defective. Contact Service as soon as possible.

Attention: No protective glass detected in the device. The device may be damaged if it is activated without the protective glass.

Symptom	Possible cause	Remedy
A weak or non-existent lithroscopy effect or effect in the tissue can be identified.	There are rings on the beam emitted by the laser fiber (Fig. 61).	Replace the protective glass as described in Section 7.1.

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8.3 Technical data

Laser type	pulsed Ho:YAG laser
Wavelength	2080 nm
Output	P ₀ = max. 35 Watt, P _p = max. 10 kW
Repetition rate	4/6/8/10/12/15/17/20/22/ 25/27/30 Hz
Pulse duration	90 ... 2000 µs
Energy	0.2/0.4/0.6/0.8/1.0/1.2/1.4/ 1.6/1.8/2.0/2.5/3.0/3.5/ 4.0 J (max. 4.8 J)
Fiber types	Bare fiber, 230 µm, 365 µm, 600 µm fiber
Fiber connection	Modified SMA fiber connector
Divergence	0.22 mrad half-angle
Aiming beam	< 5 mW, 532 nm
Power supply	Universal power supply unit 100 ... 240 VAC
Power frequency	50/60 Hz
Current consumption	max. 15 A
Cooling	Integrated water cooling system with water/air heat exchanger
Dimensions	610x300x430 mm
Weight	35 kg (incl. 2kg cooling water)
Storage/transport conditions	
Temperature	-10 °C ... 60 °C
Humidity	5% ... 95% (rel. humidity, non-condensing)
Operating conditions	
Temperature	18 °C...28 °C
Humidity	30% ... 70% (non-condensing)
Atmosph. pressure	700 ... 1060 hPa

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**8.3.1 Standard compliance
(for UL 300)**

According to IEC 60601-1, IEC 60601-2-22:

- Type of protection against electric shock:
Protection Class I
- Degree of protection against electric shock:
Applied part type B ⬆

According to IEC 60825-1:

- Laser class: 4

According to IEC 60601-1-2:

Please observe the information on electromagnetic compatibility in Section 10 (p. 65-78).

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8.4 Technical documentation

On request, the manufacturer will provide available parts lists, descriptions, instructions for adjustments, and other documentation deemed useful for suitably qualified user personnel authorized by the manufacturer when repairing equipment components that have been designated as repairable by the manufacturer.

Supply of such technical documentation relating to the device shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the device.

The maintenance and repair operations described in this manual are exempt from this rule.

① **NOTE:** We reserve the right to make design modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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9 Replacement parts/ accessories

Item	Cat. no.
Laser safety goggles Ho:YAG Laser 2080, for users with and without glasses, absorbing plastic filter, flexible temples, lightweight	27750095
Safety goggles surgical LASER, sporty version for users who don't wear glasses	UL021
Reusable fiber sets CALCULASE III laser fiber 230 µm, reusable, sterile, pack of 6 CALCULASE III laser fiber 365 µm, reusable, sterile, pack of 6 CALCULASE III laser fiber 600 µm, reusable, sterile, pack of 6	UL007 S5 UL008 S5 UL009 S5
Fiber sets for single use CALCULASE III laser fiber 230 µm, single use, sterile, pack of 6 CALCULASE III laser fiber 365 µm, single use, sterile, pack of 6 CALCULASE III laser fiber 600 µm, single use, sterile, pack of 6	UL004 S5 UL005 S5 UL006 S5
CALCULASE III set of keys	27750235
CALCULASE III interlock bridge	UL101
Fiber stripper set	27750280
CALCULASE III footswitch	27750124
Accessories for filling and draining the cooling system	UL102
Laser fiber wire tray	39555 B
Equipment trolley for the CALCULASE III	UG210

CAUTION: Never operate with defective accessories or a damaged power cable.

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Electromagnetic Compatibility (EMC) information

10 Electromagnetic Compatibility (EMC) information

The CALCULASE III UL300 is suitable for use in professional healthcare facility environment.

Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

! WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

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Accessories and cables used for EMC Compliance

Cable type	Shielded	Length [m]	Ferrite	Used for
Footswitch with cable	Yes	3	No	Activating the unit functions
Power cable 20A with primary power plug	No	5	No	Power supply

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! WARNING: The use of an accessory or cable with the CALCULASE III other than those specified in this manual may result in increased emissions or decreased immunity of the CALCULASE III. When using an accessory or cable with the CALCULASE III other than those specified in this manual, it becomes the responsibility of the user of the CALCULASE III to determine compliance with IEC 60601-1-2 when using this item.

! WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the CALCULASE III including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

i NOTE: The emission characteristics of this equipment make it suitable for use in professional healthcare environment as well as residential environment (CISPR 11 Class B). This equipment offers adequate protection to radio communication service. In the rare event of interference to the radio communication service, the user might need to take mitigation measures, such as relocating or re-orienting equipment.

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Table 1 – Compliance level for Immunity tests

Guidance and manufacturer's declaration – electromagnetic immunity

The CALCULASE III is intended for use in the electromagnetic environment specified below. The user of CALCULASE III should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV Contact discharge ± 15 kV Air discharge	± 8 kV Contact discharge ± 15 kV Air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1 kV ± 1 kV for input and output lines 100 kHz repetition	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage Dip:</u> Dip to 0 % for 1 cycle @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	<u>Voltage Dip:</u> Dip to 0 % for 1 cycle @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the CALCULASE III requires continued operation during power mains interruptions, it is recommended that the CALCULASE III be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50 Hz/60 Hz	30 A/m at 50 Hz/60 Hz	If image distortion occurs, it may be necessary to position the CALCULASE III further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz *See table 2 for wireless Proximity RF field test levels	3 V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3 V _{rms} on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V _{rms} in ISM bands	3 V _{rms} on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V _{rms} in ISM bands	

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Table 3 – Test levels for radiated and conducted immunity tests
Guidance and manufacturer's declaration – electromagnetic immunity

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Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the CALCULASE III, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b .
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz Interference may occur in the vicinity of equipment marked with the following symbol: 

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects, and persons.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CALCULASE III is used exceeds the applicable RF compliance level above, the CALCULASE III should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the CALCULASE III.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 – Emission class and group Guidelines and manufacturer's declaration – electromagnetic emissions		
The CALCULASE III is intended for use in the electromagnetic environment specified below. The customer or user of the CALCULASE III should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The CALCULASE III uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions as per CISPR 11	Class B	The CALCULASE III is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	N/A	
Voltage fluctuations/flicker emissions as per IEC 61000-3-3	Complies	

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Table 5

Recommended separation distances between portable and mobile RF communications equipment and the CALCULASE III

The CALCULASE III is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the CALCULASE III can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CALCULASE III as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@kseal.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Carneros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kseal.com

KARL STORZ Endoscopy México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ Tra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorz@kscuba.com

KARL STORZ Endoscopy México S.A. de C.V.
Av. Ejercito Nacional No. 453 Piso 2,
Colonia Granada, Alcaldia Miguel Hidalgo
C.P. 11520 Ciudad de México
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopy Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA – Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoscopy Norge AS
Stamveien 1
1483 Hagan, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
127 39 Skärholmen, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se
KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestučio st. 59 / Lenkioji st. 27
08124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 67 000
E-Mail: info-lt-n@karlstorz.com
KARL STORZ Endoscopy Danmark A/S
Skovlytten 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopy Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopy France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.com

KARL STORZ Endoscopy Iberica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopy Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Croatia d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Cesta v Gorice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Polska Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

KARL STORZ Endoszkóp Magyarország Kft.
Toberek utca 2. sz. 17/b
HU-1112 Budapest, Hungary
Phone: +36 195 096 31, Fax: +36 195 096 31
E-Mail: info-hu@karlstorz.com
KARL STORZ Endoscopy Romania srl
Str. Prof. Dr. Anton Colorian, nr. 74, Sector 4
041393 Bukarest, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
E-Mail: info-ro@karlstorz.com

KARL STORZ Endoscopy Greece M.E.P.E.*
Patriarhou Grigoriou E' 34
54248 Thessaloniki, Greece
Phone: +30 2310 304868, Fax: +30 2310 304862
E-Mail: info-gr@karlstorz.com

*Repair & Service Subsidiary

KARL STORZ Industrial**
Gedik Is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No: 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030

**Sales for Industrial Endoscopy

000 KARL STORZ Endoscopy – WOSTOK
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: info-ru@karlstorz.com

TOV LLC KARL STORZ Ukraine
Avenue Geroyiv Stalingrada Str. 2D, office 717
Kyiv, 04210/Ukraine
Phone: +38 095 000-895-0, +38-097-000-895-0,
+38 073 000-895-0
E-Mail: marketing@karlstorz.com.ua

KARL STORZ SE & Co. KG
Representation Office
Sabir Orudschow 1184, apt. 23
1025 Baku, Azerbaijan
Phone: +99 450 613 30 60
E-Mail: info-az@karlstorz.com

KARL STORZ ENDOSKOPE –
East Mediterranean and Gulf (Offshore) S.A.L.
Spark Tower 1st floor
Charles Helou St., Horch Tabet – Sin El Fil
Beirut, Lebanon
Phone: +961 1 501105, Fax: +961 1 501950
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
Saryarka, 6, BC "Arman", off. 910
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai – United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road,
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Interchange 21 Tower, Level 33
399 Sukhumvit Road,
North Klongtoey, Wattana,
10110 Bangkok, Thailand
Phone: +66 2 660 3669
E-Mail: info-th@karlstorz.com

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, M Plaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd.
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch
Room 701A Building 5 & Room 501 Building 7,
No. 3000 Longdong Avenue, Pilot Free Trade Zone,
201203, Shanghai, People's Republic of China
Phone: +86 21 60339888, Fax: +86 21 60339808
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
Room 803-805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenyang District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch
Room 02B & 03 & 04A, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 87321281, Fax: +86 20 87321286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Indonesia
Sinar Mas MSG Tower Level 37
Jl. Jend. Sudirman No. Kav. 21
Jakarta Selatan
DKI Jakarta 12920
E-Mail: infoindonesia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info-jp@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
68 Waterloo Road, Macquarie Park NSW 2113
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: karlstorz@karlstorz.com.au

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