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Does not apply to USA and Canada!

US and Canadian customers, please use the separate US instructions!

1 Article numbers

UL007, UL008, UL009

2 Intended use

The optical fibers (laser fibers) are designed for use in connection with the Calculase III laser units from KARL STORZ in the wavelength spectrum of 2100 nm.

They are intended for use in operating rooms and doctors' offices for human lithotripsy applications and for the cutting, ablation and vaporization of tissue.

2.1 Indications

The use of reusable laser fibers is indicated if the following is indicated in the opinion of the attending physician.

Lithotripsy:

Stones of all types in the ureteral system (urethra, ureter, bladder and kidneys), bile duct, gall bladder, salivary duct and salivary glands.

Urology:

- Holmium laser ablation – HoLAP
- Holmium laser resection of the prostate – HoLRP
- Holmium laser enucleation of the prostate – HoLEP
- Obstruction of the bladder neck
- Prostate and bladder incisions
- Hemostasis, vaporization and excision for treatment of benign prostate hypertrophy

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**Calculase III Laser Fiber
(reusable)**

- Superficial and invasive bladder, urethral and ureter tumors
- Condylomata
- Strictures and stenoses in the urinary tracts
- Ureteroscopic laser papillotomy
- Laser endoureterotomy
- Laser endopyelotomy
- Laser urethrotomy
- Urothelial tumors of the upper urinary tracts (UUT-UT)
- Treatment of ureteropelvic junction obstruction
- Multittract MPCNL
- Primary or refractory urethral stricture with holmium laser endourethrotomy (HLU)

2.2 Contraindications

The use of the laser is contraindicated if, in the opinion of the attending physician, the surgical method as such is contraindicated or if the patient is not able to undergo surgery or anesthesia due to his or her general condition.

Use is contraindicated if, in the opinion of the attending physician, the device is not compatible with successful completion of the planned intervention due to its technical design.

Lasers must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

Improper use (other than the intended use) is prohibited for safety reasons.

Unauthorized conversions or modifications to the laser fiber are not permitted for safety reasons.

3 Safety instructions

WARNING: Before use, check the expiration date and ensure that the packaging is undamaged. If the expiration date has passed or the packaging is damaged, the product must not be used.

WARNING: Risk of injury: Be sure to wear protective glasses tailored to the laser used. To this end, compare the details on the protective glasses with the details in the manual for the relevant laser unit.

WARNING: Risk of injury: Damaged laser fibers can lead to injuries to the patient. Laser fibers and accessories used in combination must be checked immediately before and after each application to make sure they are complete, free from damage and in full working order. Broken laser fiber parts must not be left inside the patient.

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**Calculase III Laser Fiber
(reusable)**

CAUTION: Risk of material damage when inserting the laser fiber into endoscopes. Only insert the laser fiber when the working channels are straight.

WARNING: Avoid extreme bending or tight winding of the fiber. This could cause the fiber to break. The minimum bending radius is between 1 and 3.5 cm, depending on the fiber (230 µm fiber: 1.5 cm, 365 µm fiber: 2.0 cm, 600 µm fiber: 3.5 cm). Ensure that this is also taken into consideration for storage.

WARNING: Risk of fire: Laser light from Class 4 lasers represents a source of ignition. When using the laser fiber, take suitable measures to prevent laser-induced fires and explosions.

WARNING: Risk of injury if the maximum output and energy of the laser fiber is exceeded. Only use laser fibers with the specified laser units.

Art. no.	E _{max} [J]	P _{max} [W]
UL007	2	20
UL008	4	35
UL009	4	35

4 Handling

It is recommended that the suitability of the devices for the intended procedure be checked prior to use.

- The original packaging must not be opened until immediately before the product is used.
- Check the packaging for damage: Damaged packaging may compromise the sterility of the product.
- When removing the product from the primary packaging, the rules of aseptic handling must be followed.
- The product must be examined for damage such as cracks, deformations, etc. Defective products must not be used under any circumstances.

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4.1 Remove the laser fiber from the packaging

1. Open the sterile packaging.
2. Under aseptic conditions, carefully remove the laser fiber together with the laser fiber carrier from the sterile packaging.
3. Remove the distal laser end from the carrier and carefully unwind the laser fiber from the carrier.
4. Remove the connector from the holder on the carrier.
5. Unscrew the cap from the SMA connector.

CAUTION: Risk of material damage: Do not touch the light entrance surface.

CAUTION: Risk of material damage: Do not damage the proximal light fiber end.

NOTE: Avoid contact with the fiber end and do not soil the fiber end.

4.2 Plug the connector into the laser unit

1. Remove the protection cap on the fiber SMA connector.
2. Insert the fiber connector axially into the fiber socket on the laser unit.
3. Tighten the connector with the nut by hand to fix in place.

CAUTION: Risk of material damage: Do not tilt the laser fiber connector. Do not use tools.

4.3 Functional test

Check optical laser fiber transmission.

1. Switch the laser unit on.
2. Connect the laser fiber.

The aiming beam is now available.

Calculase III = green.

3. Check the distal reproduction of the aiming beam light.

Optical laser fiber transmission is OK

Optical laser fiber transmission is **OK** if the aiming beam produces a round, evenly bright mark at a distance of 5 to 10 cm from the distal fiber end.

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Calculase III Laser Fiber (reusable)

Optical laser fiber transmission is not OK

Optical laser fiber transmission is **not OK** if the mark produced is not round, is unevenly bright or is extremely irregular (comet tail).

1. In this case, check the proximal laser fiber surface and trim the laser fiber.

If you are still unable to achieve a round aiming beam, the laser fiber must be replaced if treatment results prove inadequate.

4.4 Remove the laser fiber

1. Switch the device to standby mode.
2. Unscrew the laser fiber from the laser system. Screw the protection cap onto the fiber's SMA connector immediately afterwards.

5 Reprocessing

5.1 General warnings

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

CAUTION: When preparing and using solutions, follow the chemical manufacturer's instructions, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: Danger of damage to medical devices: The use of chemicals that have not been approved by KARL STORZ or chemical manufacturers may cause damage to the medical devices. Only use chemicals approved for the medical device for reprocessing. A helpful guide is available on the Internet at www.karlstorz.com.

CAUTION: National laws and regulations must be observed.

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Calculase III Laser Fiber (reusable)

CAUTION: Risk of material damage: Make sure the fiber connector is always sealed using the protection cap and does not touch or damage the fiber end.

NOTE: The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' (Art. no. 96216003) can be downloaded or requested by visiting www.karlstorz.com.

5.2 Accessories

Accessories required for carrying out reprocessing:
Brushes: 27652

5.3 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

5.4 Manual precleaning Brushing the surfaces

Any visible contamination or heavy soiling must be removed by cleaning the surfaces under cold running water with the aid of a brush or sponge.

Auxiliary pretreatment with ultrasound:

To aid preliminary cleaning, an ultrasound treatment must be performed. The exposure time must be at least 10 minutes at a frequency of approx. 35 kHz. The device must then be rinsed with cold water to ensure neutralization.

5.5 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual, 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003):

Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A_0 value must be taken into account when using this method.

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Calculase III Laser Fiber (reusable)

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

i NOTE: If necessary, the instrument must be dried off manually afterwards.

5.6 Assembly, inspection and care

1. The protection cap must be screwed on immediately after unscrewing the fiber from the laser (the fiber end must not be touched in the process).
 2. The protection cap must be screwed on during cleaning, and especially for ultrasound cleaning.
- The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:
- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
 - Damaged or corroded medical devices must be withdrawn from use.
 - A functional check must be carried out before using the fiber.

i NOTE: During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

5.7 Packaging systems

Only standardized and approved packaging materials and systems may be used (EN 868 Parts 2 – 10, EN ISO 11607 Parts 1 + 2, DIN 58953).

5.8 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003). The method must be selected taking into account the respective applicable national requirements and in consultation with the device and product manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

Steam sterilization using the fractionated prevacuum procedure

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i NOTE: The protection cap must be unscrewed prior to sterilization.

5.9 Limits of reprocessing

The end of the device's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use. The device's service life is limited to a maximum of 9 reprocessing cycles.

6 Description of the fiber

Depending on the fiber type, various pulse energies and medium laser outputs can be applied (see instruction manual of the Calculase III laser units).

6.1 Trimming of the fiber

The fibers are quartz fibers with a protective plastic exterior. They can be trimmed intraoperatively.

1. Remove the protective jacket using the respective fiber stripper from the fiber trimming set.
2. Trim the tip of the fiber using the ceramic knife (see Calculase III instruction manual).

! **CAUTION:** Risk of material damage: Use the sterilizable KARL STORZ fiber trimming set when handling laser fibers for intraoperative reprocessing of the distal end of the laser fiber.

i NOTE: Stripping length 10 – 15 mm.

7 Application/storage conditions

The device must be stored in a clean and dry place away from direct sunlight.

	Temperature	Rel. humidity
Storage	-20 °C...+60 °C	10 %...90 %
Operation	+10 °C...+40 °C	30 %...70 %

8 Directive compliance

This medical device bears a CE mark in accordance with the Medical Devices Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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REF UL007, UL008, UL009

2 Intended use

The optical fibers (laser fibers) are designed for use in connection with the Calculase III laser units from KARL STORZ in the wavelength spectrum of 2100 nm. They are intended for use in operating rooms and doctors' offices for human lithotripsy applications and for the cutting, ablation and vaporization of tissue.

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
Urology:

- Holmium laser ablation – HoLAP,
- Holmium laser resection of the prostate – HoLRP,
- Holmium laser enucleation of the prostate – HoLEP,
- Obstruction of the bladder neck
- Prostate and bladder incisions
- Hemostasis, vaporization and excision for treatment of benign prostate hypertrophy
- Superficial and invasive bladder, urethral and ureter tumors

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Calculase III Laser Fiber (reusable)

- Condylomata
- Strictures and stenoses in the urinary tracts
- Ureteroscopic laser papillotomy
- Laser endoureterotomy
- Laser endopyelotomy
- Laser urethrotomy
- Urothelial tumors of the upper urinary tracts (UUT-UT)
- Treatment of ureteropelvic junction obstruction
- Multitract MPCNL
- Primary or refractory urethral stricture with holmium laser endourethrotomy (HLU)

2.2 Contraindications

The use of the laser is contraindicated if, in the opinion of the attending physician, the surgical method as such is contraindicated or if the patient is not able to undergo surgery or anesthesia due to his or her general condition. Use is contraindicated if, in the opinion of the attending physician, the device is not compatible with successful completion of the planned intervention due to its technical design.

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NOTE: If necessary, the instrument must be dried off manually afterwards.

5.6 Assembly, inspection and care

1. The protection cap must be screwed on immediately after unscrewing the fiber from the laser (the fiber end must not be touched in the process).
2. The protection cap must be screwed on during cleaning, and especially for ultrasound cleaning.

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- A functional check must be carried out before using the fiber.

NOTE: During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

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The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) at 132 °C – 137 °C with a minimum exposure time of 3 minutes and a maximum of 18 minutes.

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NOTE: The protection cap must be unscrewed prior to sterilization.

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1. Remove the protective jacket using the respective fiber stripper from the fiber trimming set.
2. Trim the tip of the fiber using the ceramic knife (see Calculase III instruction manuals).

CAUTION: Risk of material damage: Use the sterilizable KARL STORZ fiber trimming set when handling laser fibers for intraoperative reprocessing of the distal end of the laser fiber.

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