
Surgical Dr.Lite 2066 (2066-05, 2066-06, 2066-07)

Installation/Operating Instruction

 **BenQ Medical Technology Corporation**
7th Floor, No. 46, Zhou-Z Street, Nei-Hu,
Taipei 114, Taiwan
Customer Service e-mail :
Service@BenQMedicalTech.com
BenQMedicalTech.com
Customer Service Hotline :
Taiwan : 0800-015-533
China : +86 21-6327-7161~3 Ext. 812
International : +886 2-2601-7227 Ext. 5932

This operating instruction is intended to provide users with the proper information required for correct operation and maintenance of BenQ surgical light. In order to avoid danger to patient or users and/or damage to the equipment, please review the contents of this instruction thoroughly and carefully before use.

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BenQ STANDARD LIMITED WARRANTY

BenQ warrants all new surgical light sold by it directly or through a dealer or other authorized representatives, with expectation to replacement parts, spares, pads, and accessories, to be free from defects in material or workmanship, under normal use and service, for a period of one year. The warranty period shall begin with one year after shipment from BenQ.

BenQ's responsibility and liability shall be limited to the replacement of any parts which we, BenQ, determine to be defective within the applicable warranty period, labor will not be included in the warranty coverage. No dealer, agent, or employee is authorized to extend or enlarge this warranty.

EC REP

ARTMED Ltd.

H-1138 Budapest, Viza u 7/c, B303, Hungary.

TABLE OF CONTENTS

/I/ Summary of Safety and General Precautions	Page 2
Section 1 : Section Introduction	
1.1 General	Page 4
1.2 Complete set	Page 5
1.3 Power Source Requirement	Page 5
1.4 Explanation of Labels	Page 6
1.5 Cleaning or Disinfecting	Page 8
Section 2 : Installation Introduction	
2.1 Ceiling Mounting Methods	Page 9
2.2 Installation Procedure of Suspension Arm and Light-head	Page 12
2.3 Power Connection Instructions for Ceiling Power Box	Page 15
2.4 To exchange fuse	Page 16
2.5 The specification of fuse	Page 16
Section 3 : Operating Introduction	
3.1 Wall Control Box (Option)	Page 17
3.2 Control Panel	Page 18
3.3 Endo. Light (optional)	Page 20
3.4 Initial Status of System Start	Page 20
3.5 Focused Pattern Size Adjustment and Light-head Positioning	Page 20
3.6 Brake/Friction Adjustment	Page 21
3.7 Spring Arm Adjustment	Page 23
3.8 Bulb Change	Page 24
Section 4 : Adjustment and Maintenance	
4.1 Troubleshooting	Page 25
Section 5 : Technical Data	
5.1 Technical Specifications	Page 26
5.2 Range of Circular Motion	Page 28
5.3 Environmental conditions for storage and transportation	Page 31

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WARNING:

PINCHING HAZARD: Pinch points are composed in extreme light-head articulations. User should be aware of and examine these points before operating the light.

- ⊙ Potential pinch and/or crush hazard exists if hands or fingers are placed between edges of yoke and light-head during extreme spring arm posturing.

FAILURE HAZARD:

- ⊙ Do not use BenQ light parts or bulbs for other than their stated function or parts or bulbs manufactured by other manufacturers to avoid possible patient injury, user injury, or surgical light or accessory failure. Any warranties could also be rendered void by the use of non-original parts or bulbs.
- ⊙ For safe and reliable operation of the lamp, installation, maintenance and repair of the equipment must be done only by qualified service personnel. Preventive maintenance of the lamp requires a regular schedule.
- ⊙ Serious injury or lamp damage may result if user of the lamp is not fully familiar with the controls for lamp operation and needed positioning.

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EXPLOSION HAZARD: Do not use in the presence of flammable anesthetics.

BURNT HAZARD:

- ⊙ Do not change bulb when the cover of light-head is still warm. Quartz Halogen Bulbs generate considerable heat, which is emitted through the top of the light-head. The covers will be quite warm to be touched.

PERSONAL INJURY HAZARD:

- ⊙ Do not remove or modify light fixture by unauthorized technicians. Refer servicing to authorized technicians.
- ⊙ Do not attempt to adjust light beam focus or suspension arm assembly. Refer servicing to authorized technicians.
- ⊙ Use a lift capable to safely lower or install the suspension arm assembly, spring arm, and light-head/yoke assembly.
- ⊙ Make sure the optical core handle assembly of light-head is supported before operate the light.

- ⊙ Do not work under unsupported light fixture or load.
- ⊙ Do not operate the light if the Heatsorb filters or cover lens are broken or removed. Operation without the Heatsorb filters or cover lens or with broken Heatsorb filters or cover lens could result in injury to the patient or user.

⚠ INTERFERENCE HAZARD:

- ⊙ Do not use surgical light in the presence of magnetic field or nearby other interference equipment.

⚠ DISPOSAL OF HAZARDOUS WASTE EQUIPMENT OR SUBSTANCE:

- ⊙ You must dispose of your waste equipment or substances by handing it over to or collected by the licensed company or organization which is eligible to handle with hazardous material. Proper disposal of materials will help to conserve natural resources and ensure that it is recycled or handled in a manner that protects human health and environment.


⚠ CAUTION:

EQUIPMENT DAMAGE:

1. Do not bump light-head or suspension arm assembly into another light-head, walls, or other equipment.
2. Prevent leakage of fluids into inside of light-head.
3. Do not expose the lamp to excessive moisture.
4. Do not use harsh cleaning, disinfecting, or anti-static liquid on the light. Use of phenolic-, iodophor-, glutaraldehyde-, alcohol-based disinfectants on surface of the light could cause some degree of staining, pitting, and/or discoloration damage on painting or optical cover.
5. Use only clean, soft, and lint-free cloth to wipe both internal or external surfaces of this light when cleaning.
6. Do not touch Quartz Halogen Bulb with bare hands to prevent damage to the lamps and decrease to the life of the bulbs.
7. Do not touch optical reflector and Heatsorb filter when changing bulbs.
8. The ceiling mounting disc plate must be level in two directions. Improper operation of the suspension arm assembly could result in equipment damage.

SECTION 1 : SECTION INTRODUCTION

1.1 GENERAL

BenQ's Dr.Lite 2066 surgical lights series are a single evolutionary reflector light with numerous facets to provide a field of illumination in surgical, diagnosis, or treatment applications. The single light-head model is intended for general and large surgical, diagnosis, and/or treatment applications and may be used in the situation, where the optimal illumination, cold light, color temperature and color rendering index suitable for procedures observations can be provided where blood vessels, organ tissues, and nerves of surgical site can be differentiated. The color temperature of Dr.Lite is 4200 ± 200 °K that is close to average day light, operating rooms for all surgical procedures and where a safe condition is required.

This lighting condition will come in handy to all surgeons as if all the surgical procedures are performed at 8 o'clock in the morning of summertime. Dr.Lite 2066 surgical lights series deliver cold light required for surgical procedures where blood coagulation and tissue dryness must be prevented and provide surgical team with good working environment and eliminate chances of uncomfortable illumination owing to overheat. Automatic back-up bulb system will detect a burnt bulb and activate alternative working bulb to ensure supply of illumination throughout entire surgical procedure. The operation life of the surgical light under normal use is 10 years.

1.2 CONTENTS OF WOODEN CASE

1. Suspension Arm Assembly - not more than 3 pcs.;
2. Spring Arm - not more than 3 pcs.;
3. Mounting plate - 1 pcs.;
4. Ceiling housing - 1 pcs.;
5. Mounting kit - 1 pcs.;
6. Power Box - 1 pcs.;
7. Light-heads - 1 pcs.;
8. Lighting unit with built-in CCD - 1 pcs.;
9. Optical Core Handle Assembly (with Sterilizable Handle) - 1 pcs.;
10. Remote control - 1 pcs.;
11. Separate CCD - 1 pcs.;
12. Bracket for mounting a camera - not more than 1 pcs.;
13. TFT-display - not more than 3 pcs.;
14. Control unit wall - 1 pcs.;
15. Medical PC ORLIS - 1 pcs.;
16. Handle with built-in LEDs EndoLED - 1 pcs.;
17. Vertical Gimbals Joint - not more than 3 pcs.;
18. Halogen lamp - not more than 10 pcs.;
19. LEDs - not more than 160 pcs.;
20. AAA batteries - not more than 2 pcs.;
21. Power cable - not more than 5 pcs.;
22. Instruction for use - 1 pcs.;
23. Cover Lens - not more than 1 pcs.;
24. Extension pipe - not more than 2 pcs.

1.3 POWER SOURCE REQUIREMENT

1. Dr.Lite series surgical lights can accept applied power source 100-240Va.c, 50/60Hz. (Please see the label on the power box before connecting power source.)
2. If there is a need to extend the length of power cable between the power box and surgical light, please follow the following formula to calculate the precise cross section area for each individual core.

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A= The cross section area of the extended power cable (mm²).

L= The length of extended power cable between power box and surgical light (m).

P= Power consumption (W), please refer the power consumption to the light-head label.

CAUTION:


Always limit the applied power source at the range of $\pm 5\%$ to avoid equipment damage.

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
【 LABEL ATTACHED PLACE AND EXPLANATION 】

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1.5 CLEANING OR DISINFECTING

1. To remove spots or stains, wipe off with a lightly dampened cloth and ordinary disinfectants after using Dr.LiTE 2066 Surgical Light Series.
2. To sterilize maneuvering handle of light-head, the recommended parameters through steam sterilization are:

Steam sterilization at 121°C ; 1.1bar ; 15 minutes.

Steam sterilization at 134°C ; 2.1bar ; 5 minutes.

CAUTION:

Use only clean, soft, and lint-free cloth to wipe both internal or external surfaces of this light when cleaning.

SECTION 2 : INSTALLATION INTRODUCTION

2.1 CEILING MOUNTING METHODS

✓ **Figure(2.1.1):** The distance between the lowest point of spring arm clevis and the finished floor is at least 6'6"(200cm) after completion of installation. Our Dr. Lite lights provide 3 different floor-height for you to choose, such as 8'10"(270cm), 9'4"(285cm), 9'10"(300cm).


Figure (2.1.1)

Situation (1) When the distance existing ceiling between structure ceiling and floor is within 9'10" (300cm), the installation of mounting disc plate onto existing ceiling structure is shown as figure (2.1.2).

Studs : 3/4"-10UNC thread(or M18), 4 pieces.



Figure (2.1.2)

Mounting Disc Plate

- ✓ If a ceiling height exceeds 9'10" (300cm), a properly reinforced intermedium frame is required to support the lamp. If there is a finished ceiling, a properly reinforced intermedium frame is also needed for installation of the lamp. The following is the methods of a properly reinforced intermedium frame, depending on ceiling height.

Situation (2) If the distance between existing ceiling structure and finished ceiling is within 11.8" (30cm), the installation instruction of a properly reinforced intermedium frame shown as figure (2.1.3). Use expansion (foundation) bolts to fix the 4 holes, $\phi 0.87"$ (22mm) of mounting disc plate in place, then use 4 connect nuts and studs to extend your needed length. Ceiling plate must be level in two directions. The upper surface of the mounting disc plate must align with the finished floor shown as figure (2.1.3).


Figure (2.1.3)

- 1: Put M16 expansion (foundation) bolts into existing ceiling structure for 3.15" (8cm).
- 2: 3/4"-10UNC studs (or M18) studs.

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
Figure (2.1.4)

Reinforced tubes:
 $\phi 1.22"$ (31mm, inner) $\times \phi 1.5"$ (38mm, outer).

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
Figure (2.1.5)

✓ **Manufacturing of large-tube extension frame:**

- Length of large tube between each point on both cutting edges must be the same. This will ensure that frame is always align and level after weld with the steel plates on each end shown as figure (2.1.6).








Figure (2.1.6)

Caution : The maximum length of extension tube is no longer than 800mm.


Improper electrical grounding or no grounding between extension tube and disk may result in personal injury or equipment damage.

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1  Mounting Disc Plate Locking Bolt Suspension Arm Assembly Fix suspension arm assembly onto mounting disc plate with 6 locking bolts. Note: Always level mounting disc plate before fixing suspension arm assembly.	2  Remove the Set Screw on the Mount of Spring Arm.
3  Insert the Male Connector into Mount on Spring Arm, and fixed the Male Connector with Set Screw.	4  Insert the Mount of Spring Arm onto Extension Arm, and fixed with Set Screw.
5  Align braking setscrew hole with a concavity on inner plastic tube precisely and screw a braking screw properly to have proper brake on related spring arm.	6  Use screwdriver to properly tighten braking setscrew, and do micro-adjustment after light-head is installed.

• Remove screw on spring arm with phillips screw driver.	• Lift safety sleeve slightly, twist to the right, and then pull right up. • Remove safety segment.
• Shown as the figure, Separating connector and light head.	• Shown as the figure, combine with connector and spring arm.
• Insert spigot of light-head yoke into socket on spring arm and remove self-adhesive label.	• Replace safety segment back to socket. • Make sure that safety segment has been inserted into the socket.

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

14

- Attaching this label to spring arm sleeve.
P.S : Only qualified technician or engineer should be allowed to peel off the label once label has been applied.


15

- Insert optical core assembly into light-head.
- Pay attention upon 4 electrical sockets when connecting.


16

- Twist slightly and place ceiling shroud to cover central base of suspension arm assembly.

17

- Twist and turn the ceiling shroud again to pass upper suspension arm.

18

- Lock ceiling shroud with setscrews.

2.3 POWER CONNECTION INSTRUCTIONS FOR POWER BOX

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Installation Procedures:

1. Place and fix the power supply box in the ceiling mounting frame (top of mounting disc plate) or other proper position. Shown as figure (2.3.1), connect power cord (C) to applied power 100-240 Va.c.. Turn on power switch (A) and check if green light is lit. If yes, the power supply is normal. Please turn off power supply and do next step. If no, check power supply.
2. Shown as figure, connect Vd.c. cable with power cable of surgical light. Turn on power switch on the box and test the surgical light system. If the intensity LED is at LED 3, this means that the light system is working normally. If no, please check the power cable and do trouble shooting process.
» Male connector and female connector 4P is connected for RS-485 system (Wall control box system, ORLIS).
3. After complete power connection and installation, turn on power switch on the box and disconnect the applied power supply to wire all cables and connectors properly.

(I) Suspension Arm
Identification Label-1

(J) Suspension Arm
Identification Label-2

(K) Suspension Arm
Identification Label-3


The mark of primary light head is 『Lamp (1)』.
The mark of secondary light head is 『Lamp (2)』.
The mark of third light head is 『Lamp (3)』.
The mark of the carry arm of TFT is 『TFT』.
The mark of AC power in is 『AC 100- 240Va.c.』.

WARNING:

- Improper electrical grounding or no grounding may result in personal injury or equipment damage.
- To check the mark, before connecting with Vd.c. and power cable.
- Each light head has its own power supply box. This means that dual light head surgical light has 2 power supply boxes.


2.4 To exchange fuse


Use screw driver to loose the Fuse cover.

Replace the burn out Fuse with the same specific Fuse.


Reinstall the Fuse cover.


Use screw driver to lock the Fuse cover and then finish the Procedure.

2.5 The specification of Power Supply and Fuse


SECTION 3 : OPERATING INTRODUCTION


3.1 Wall control box system (Option)

Single Lamp system

Dual Lamp system

Triple Lamp system

3.2 CONTROL PANEL


Dr.LiTE 2066 series is equipped with light-head mounted control panel to control system on/off and intensity adjustment and to show intensity level and bulb status.

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
© Dr.LiTE 2066 light-head is mounted with a film-keypad control panel for:

【Control Panel】	【CCD Control Panel (optional)】
(A) Power ON/OFF. (B) The power "ON" Indicator. (C) Intensity of Illumination. (Increase/Decrease) (D) Burnt bulb monitoring for primary and secondary bulb. (E) Indicators for intensity level. (F) Halogen bulb control selector. (G) Endo. Light LED control selector. (H) Penetrating Light LED control selector.	F. CCD ON/OFF. G. The CCD "ON" Indicator. H. Auto/Manual Focus control switch. I. Auto Focus "ON" Indicator. J. Zoom Ratio Adjustment. (Increase/Decrease). K. Manual Focus Adjustment. (Increase/Decrease)


Control Panel


CCD Control Panel (optional)

Penetrating / Endo. Light Control Panel (Optional) :

Remote Control

- ⊙ The user presses a button (L1 、 L2 、 or L3) to select the desired lighthouse ; thus, the related LED light up to indicate the desired lighthouse on the remote control.
- ⊙ When pressing the power switch (T) to light up or light out the desired lighthouse, the remote control can't be operated within 2 seconds of power-on or power-off valid time.
- ⊙ This remote control operates with power save mode, when the user stop operation for 10 seconds, CPU will enter into dormant mode until the user presses any functional buttons to awaken CPU out of dormant mode and stay at the selected operation status prior to entering into dormant mode.
- ⊙ The remote control operation can be activated from a distance of up to 5M radius.
- ⊙ The remote control uses infrared transmission to control the surgical light. Please avoid the presence of obstacles in the emission route of remote control for interfering.


- ⊙ The infrared transmission of the remote control is directional, please target the remote control to the direction of signal receiver Sensor.
- ⊙ Use two AAA batteries, make sure to take out batteries if the batteries are not used for a long time to avoid battery-liquid leakage and erosion on circuit board.
- ⊙ Battery life of this remote control is estimated to be three months of usage, make sure to replace battery periodically to ensure normal operation.
- ⊙ Please check the battery if the remote control does not work well.
- ⊙ Push the function selector (W) 、 (V) 、 (S) to control the intensity of halogen bulb (V) 、 Endo. Light LED (W) .

3.3 Endo. Light (optional)

Dr. lite 2066 has the fade lighting for MIS surgery, and to provide three LED lighting source in sterilizing handle, as shown in figure.

The Endo. Light function be supported at Penetrating / Endo. Light control panel and remote controller.

Push the function selector (M) / (W) to star the Endo. Light LED, and using the function key (C) / (Q) and (R) to control the intensity of LED.


3.4 INITIAL STATUS OF SYSTEM START

Dr.LiTE2066 series has the "automatic system self-turn-on design" whenever applied power is connected to the surgical light system to ensure turning the system itself on automatically as applied power coming back after temporary power shortage.



3.5 FOCUSED PATTERN SIZE ADJUSTMENT AND LIGHT-HEAD POSITIONING




- ⊙ By turning the sterilizable maneuvering handle to adjust size of the focused pattern or using optional focus control panel or remote control to adjust the focused pattern size.
- ⊙ Users may reposition light-head through any positions with the sterilizable maneuvering handle or auxiliary handles located at either side of the control box.

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PROBLEM	SOLUTION
Inner hub of suspension arm assembly drifts horizontally.	Shown as right figure, use Allen wrench M6 to tighten the setscrews, by turning clockwise. 
Outer hub (spring arm) of suspension arm assembly drifts horizontally.	Shown as right figure, use screwdriver to tighten the setscrews, by turning clockwise. 
Spring Arm drifts vertically.	1. Shown as right figure, insert and push Allen wrench M5 fully into the opening to adjust weight capacity. 2. To increase weight capacity, by turning Allen wrench toward counterclockwise. 3. To decrease weight capacity, by turning Allen wrench toward clockwise. 

PROBLEM	SOLUTION
Light-head movement height is too high or too low.	1. Shown as right figure, insert and push Allen wrench M5 fully into the opening to adjust spring arm height. Our factory set-up is 45°. 2. To have higher, fixed lowest position, by turning Allen wrench toward counter- clockwise. 3. To have lower, fixed lowest position, by turning Allen wrench toward clockwise. 
Yoke contains rotation drift.	Shown as right figure, use bent nose plier to correct rotational drifting by turning the locknut inside the hub clockwise. 
Light-head contains rotational drift.	Shown as right figure, use bent nose plier to correct rotational drifting by turning the locknut inside the hub clockwise. 

※If you still can not solve problem, please contact our local agent or our service division in Taiwan with detailed problem, condition, and series number of a problem unit.

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1. Insert and push Allen wrench M5 fully into the opening **X** shown as figure(a).
2. To have higher fixed lowest position, by turning Allen wrench toward counter-clockwise.
3. To have lower fixed lowest position, by turning Allen wrench toward clockwise.

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1. Insert and push Allen wrench M5 fully into the opening **Y** shown as figure(a).
2. To increase weight bearing capacity, by turning Allen wrench toward counterclockwise.
3. To decrease weight bearing capacity, by turning Allen wrench toward clockwise.

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


Figure (a) Spring Arm Adjustment

3.8 BULB CHANGE

- Release three locking bolts to remove optical core handle stem out of light-head.


- Remove optical core handle stem out of light-head and inspect the damage condition of the primary (main) and secondary (sub.) bulb.
- Wear cotton gloves when replacing the burnt primary (main) bulb. And inspect burnt primary (main) bulb, secondary (sub.) bulb, and bulb socket for any abnormality. If the spare bulb is abnormal, replace it accordingly.

Warning:

1. Do not change bulb when the bulb of the surgical light is still warm.
2. Wear gloves to maintain the tidiness of bulb socket, please take hold of protection film, once the replacement is completed, take out protection film and use wiping cotton towel to clean bulb and bulb socket.
3. Pay attention to the position of primary/secondary bulb and replace accordingly.

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4.1 TROUBLESHOOTING

1. Use operating instructions to verify trouble symptoms.
2. After symptom has been verified, refer to the following table and select example that is most appropriate to your problem. Follow recommended correction.
3. Use electrical schematics as guide in understanding operation of the light and in locating problems.
4. Refer to Exploded Views and Parts Lists for guidance when attempting repairs.

© New Installation:

PROBLEM	SOLUTION
Bulb will not operate.	1. Use multiple-meter to check the applied power for power supply box or wall control box. 2. Check polarity of wiring between power supply box and light. Follow "Installation/Operating Instruction" to rewire or reverse polarity of wires. 3. Make sure fuse(s) is fine. 4. Check for 24 Vd.c. output at power supply box output terminal. Follow "Installation/Operating Instruction" to adjust Vd.c. to reach 24 Vd.c.. 5. Check for 22.5 -22.8 Vd.c. input at Main, Sub, and Com points on Lamp Power Board (round PCB) inside light-head. 6. Check for proper connection and normal function between light-head control panel and control board (rectangular PCB). 7. Make sure light bulb(s) is properly installed and power cable connection between bulb socket(s) and power PCB is normal.

© After Start-Up:

PROBLEM	SOLUTION
Bulb will not operate.	1. Check all items under "New Installation". 2. Check STAND BY and/or bulb failure indicator on control panel. If the indicator is not lit, power supply is abnormal. Replace new bulb(s).
Light flickers.	1. Check for loose wires in light-head, power supply box, or wall control box. 2. Check for any arcing on brushes or brush holder. 3. Defective bulb socket/holder assembly.
Low light intensity, or beam pattern, poor.	1. Bulb voltage low. (22.5-22.8 Vd.c. input at Main, Sub, and Com points on Lamp Power Board (round PCB) inside light-head.) 2. Bulb defective (glass envelope clouded, etc).
Short bulb life.	1. Bulb voltage is too high. Follow "Installation/Operating Instruction" to adjust the High Voltage down to 22.5 -22.8 Vd.c. (Volts DC.) at bulb end.
The remote control doesn't function.	1. Please check the condition of storage battery. 2. Please check if the channels (Main, Sub, Com) of the target light-head. 3. Please avoid the presence of obstacles in the emission route of remote control to interfering.

SECTION 5 : TECHNICAL DATA

5.1 TECHNICAL SPECIFICATIONS

LIGHTING DATA		2066-05	2066-06	2066-07
Dia. of Reflector		42 cm (16.5")	50 cm (19.7")	62 cm (24.4")
Max. Intensity at 1m (39.4") Distance		Max. 100,000Lux	Max. 140,000Lux	Max. 160,000Lux
Color Temperature		4200±200 °K		
Color Rendering Index		94.89		
Total Irradiance		370W/m ² (3.7mW/m ² lx)	518W/m ² (3.7mW/m ² lx)	592W/m ² (3.7mW/m ² lx)
Intensity Level Indication		8 Steps		
Focused Pattern Adjustment	Minimum	7 cm		
	Maximum	25 cm	27 cm	30 cm
Minimum Depth of Field without Refocusing		110 cm	120 cm	120 cm
MECHANICAL DATA		2066-05	2066-06	2066-07
Dia. of Light-head		50cm (19.6")	59 cm (23.2")	71.5 cm (28.14")
Min. Height from Ground to Ceiling		270cm (106.3")		
Sterilizable Temperature of Maneuvering Handle		134 °C		
Light-head Weight		10.5kg/23.0lb±5%	14kg/30.8lb±5%	16.5kg/36.3lb±5%
Total Weight	2066-05	47.0 kg±10% (103.5 lb)		
	2066-06	50.5 kg±10% (111.2 lb)		
	2066-07	54.5 kg±10% (120 lb)		
	2066-65	84.0 kg ±10% (185.0 lb)		
	2066-75	87.8 kg ±10% (193.4 lb)		
	2066-66	87.5 kg ±10% (192.7 lb)		
	2066-76	91.3 kg ±10% (201.1 lb)		
	2066-765	128.0 kg±10% (281.9 lb)		
ELECTRICAL DATA				
AC Input		100 - 240 Va.c.		
Power Consumption		Max. 150 W		
Number of Bulbs		2 (Primary and Secondary)		
Bulb		22.8 V / 150 W		
Bulb Average Life		1,000 Hours		
Applied Part			Type B, IEC 60601-1, IEC 60601-2-41	
Fuses		5Ax2, 10Ax1		
Batteries for Remote Control		AAA x 2		

【Option】 © Wall control box system


ITEM	SPECIFICATION
Wall control box Dimensions	2066-05/06/07A : 130×170×40 mm
	2066-55/65/66/75/76/77A : 260×170×40 mm
	2066-555/655/665/755/775/765/666 /777A : 390×170×40 mm
Input/output terminal	5 Pin Sub
Signal type	RS-485
Operating voltage	24Vd.c.±10%
Power consumption (working)	3.0W max.
Operating temperature and humidity	5°C - +40°C ; 20 - 80%RH
Weight(Box only)	2066-05/06/07A : 1.1 kg
	2066-55/65/66/75/76/77A : 2.2 kg
	2066-555/655/665/755/775/765/666 /777A : 3.3 kg

【Option】 © Built-in Endo. Lighting optical core assembly


Item	Specification
Max. Intensity at 1m (39.4") Distance	300Lux ±20% Lux
Color Temperature	4,300±200 Kelvin
Power Consumption	15 W
Intensity Level Indication	8 steps/electron
LED Life	20,000 hrs

5.2 RANGE OF CIRCULAR MOTION


© Dr.Lite 2066-07/06/05


© Dr.Lite 2066-77/76/75/66/65/55


© Dr.Lite 2066-666/665/655/555


© Dr.Lite 2066-07L/06L/05L

© Dr.Lite 2066-77L/76L/75L/66L/65L/55L


© Dr.Lite 2066-777L/766L/765L/666L/665L/655L/555L

5.3 Environmental conditions for storage and transportation

Operating/Shipping/Storage Environment :	
Ambient temperature	Shipping/Storage : -30~70°C
	Operating : 10~40°C
Relative Humidity	Shipping/Storage : 20~80 %
	Operating : 30~75 %
Atmospheric pressure	Shipping/Storage : 500~1060 hpa
	Operating : 700~1060 hpa

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
TriLite LED S600 Series Surgical Light

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 BenQ Medical Technology Corporation
7th Floor, No. 46, Zhou-Z Street, Nei-Hu,
Taipei 114, Taiwan

Customer Service e-mail :
Service@BenQMedicalTech.com

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H-1138 Budapest, Viza u 7/c, B303, Hungary.

TABLE OF CONTENTS

 Summary of Safety and General Precautions	_____ Page 2
Section 1 : Section Introduction	
1.1 General	_____ Page 4
1.2 Contents of Wooden Case	_____ Page 5
1.3 Power Source Requirement	_____ Page 5
1.4 Explanation of Precaution Labels	_____ Page 6
1.5 Cleaning or Disinfecting	_____ Page 8
Section 2 : Installation Introduction	
2.1 Ceiling Mounting Methods	_____ Page 9
2.2 Installation Procedure of Suspension Arm and Light-head	_____ Page 12
2.3 Power Connection Instructions for Ceiling Power Box	_____ Page 15
2.4 To replace fuse	_____ Page 16
2.5 The specification of fuse	_____ Page 16
Section 3 : Operating Introduction	
3.1 Control Panel	_____ Page 17
3.2 Endo. Light Control	_____ Page 19
3.3 ASM Control(optional only for 70 Light-head)	_____ Page 19
3.4 Initial Status of System Start	_____ Page 20
3.5 Focused Pattern Size Adjustment and Light-head Positioning	_____ Page 20
3.6 Brake/Friction Adjustment	_____ Page 20
3.7 Spring Arm Adjustment	_____ Page 22
Section 4 : Adjustment and Maintenance	
4.1 Troubleshooting	_____ Page 23
Section 5 : Technical Data	
5.1 Technical Specifications	_____ Page 26
5.2 Range of Circular Motion	_____ Page 27

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ELECTRICAL SHOCK HAZARD:

- ⊙ Do not remove cover of electrical box and/or light-head before disconnecting electrical power to the system or other than those stated in this instruction.
- ⊙ Do not attempt to clean light-head unless power is turned off and light-head is cool adequately.

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CAUTION:

EQUIPMENT DAMAGE:

- ⊙ Do not bump light-head or suspension arm assembly into another light-head, walls, or other equipment.
- ⊙ Prevent leakage of fluids into inside of light-head.

- ⊙ Do not expose the lamp to excessive moisture.
- ⊙ Do not use harsh cleaning, disinfecting, or anti-static liquid on the light. Use of phenolic-iodophor-, glutaraldehyde-, alcohol-based disinfectants on surface of the light could cause some degree of staining, pitting, and/or discoloration damage on painting or optical cover.
- ⊙ Use only clean, soft, and lint-free cloth to wipe both internal or external surfaces of this light when cleaning.
- ⊙ The ceiling mounting disc plate must be level in two directions. Improper operation of the suspension arm assembly could result in equipment damage.

Figure Illustration :

	Warning: indicating possibility of personal injury.
Model	Model
REF	Product Number
SN	Serial Number
	Grounding
	Power On/Off
	LED Failure Indicator
	Illumination Increase
	Illumination Decrease
	Endo. Light LED On/Off
	Automatic Shadow Management On/Off


SECTION 1 : SECTION INTRODUCTION

1.1 GENERAL

BenQ's TriLite LED S600 series surgical lights use LED technology to provide a field of clear illumination in surgical, diagnosis, or treatment applications. The single light-head model is intended for general and large surgical, diagnosis, and/or treatment applications and may be used in the situation, where the optimal illumination, cold light, color temperature and color rendering index suitable for procedures observations can be provided where blood vessels, organ tissues, and

nerves of surgical site can be differentiated. The color temperature of TriLite LED S600 is $4,200 \pm 200$ °K that is close to average day light, operating rooms for all surgical procedures and where a safe condition is required.

This lighting condition will come in handy to all surgeons as if all the surgical procedures are performed at 8 o'clock in the morning of summertime. TriLite LED S600 surgical lights series deliver cold light required for surgical procedures where blood coagulation and tissue dryness must be prevented and provide surgical team with good working environment and eliminate chances of uncomfortable illumination owing to overheat.

1.2 CONTENTS OF WOODEN CASE

There are cartons for suspension arm assembly and light head assembly separately inside the wooden case.

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1.3 POWER SOURCE REQUIREMENT

1. TriLite LED S600 series surgical lights can accept applied power source 110/220Vdc, 50/60Hz. (Please see the label on the power box before connecting power source.)
2. If there is a need to extend the length of power cable between the power box and surgical light, please follow the following formula to calculate the precise cross section area for each individual core.

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

1.5 CLEANING OR DISINFECTING

1. To dismount applied power before cleaning or disinfecting.
2. To remove spots or stains, wipe off with a lightly dampened cloth and ordinary disinfectants.
3. To sterilize maneuvering handle of light head, the recommended parameters through steam sterilization are:
 Steam sterilization at 121°C ; 1.1bar ; 15 minutes.
 Steam sterilization at 134°C ; 2.1bar ; 5 minutes.

⚠ CAUTION:

- Do not use harsh cleaning, disinfecting, or anti-static liquid on the light. Use of phenolic-, iodophor-, glutaraldehyde-, alcohol-based disinfectants on surface of the light could cause staining, pitting, and/or discoloration damage on painting or plastic parts.
- Use only clean, soft, and lint-free cloth to wipe both internal or external surfaces of this light when cleaning.
- Ensure that open end of sterilizable maneuvering handle is face down and the handle must not have any contacts with other objects during steam sterilization.

SECTION 2 : INSTALLATION INTRODUCTION

2.1 CEILING MOUNTING METHODS

Figure(2.1.1): The distance between the lowest point of spring arm clevis and the finished floor is at least 6'6"(200cm) after completion of installation. Our TriLite LED S600 lights provide 3 different floor-height for you to choose, such as 8'10"(270cm), 9'4"(285cm), 9'10"(300cm).


Figure (2.1.1)

Situation (1) When the distance existing ceiling Existing Ceiling between structure ceiling and floor is within 9'10" (300cm), the installation of mounting disc plate onto existing ceiling structure is shown as figure (2.1.2).

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Situation (2) If the distance between existing ceiling structure and finished ceiling is within 11.8" (30cm), the installation instruction of a properly reinforced intermedium frame shown as figure (2.1.3).

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1  Mounting Disc Plate Locking Bolt Suspension Arm Assembly • Fix suspension arm assembly onto mounting disc plate with 6 locking bolts. Note: Always level mounting disc plate before fixing suspension arm assembly.	2  • Remove the Set Screw on the Mount of Spring Arm.
3  • Insert the Male Connector into Mount on Spring Arm, and fixed the Male Connector with Set Screw.	4  • Insert the Mount of Spring Arm onto Extension Arm, and fixed with Set Screw.

 Align braking setscrew hole with a concavity on inner plastic tube precisely and screw a braking screw to have proper brake on the spring arm. Caution: Always align braking setscrew hole with a concavity on inner plastic tube precisely before screwing a braking setscrew.	 Use screwdriver to properly tighten braking setscrew, and do micro-adjustment after light-head is installed.
 Remove screw on spring arm with phillips screw driver.	 - Lift safety sleeve slightly, twist to the right, and then pull right up. - Remove safety segment.
 Insert spigot of light-head yoke into socket on spring arm.	 - Replace safety segment back to socket. - Remove protection sheet of cover lens of light-head.

 11 • Pull the safety sleeve downwards. • Screw the locking screw back to lock the sleeve with the socket.	 12 • Attaching this label to spring arm sleeve. Note : Only qualified service personnel is allowed to peel off the label once label has been applied.
 13 Ceiling Shroud • Twist slightly and place ceiling shroud to cover central base of suspension arm assembly.	 14 • Twist and turn the ceiling shroud again to pass upper suspension arm.
 15 Setscrew • Lock ceiling shroud with setscrews.	

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The mark of primary light head is 『Lamp (1)』 .
The mark of secondary light head is 『Lamp (2)』 .
The mark of third light head is 『Lamp (3)』 .
The mark of the carry arm of TFT is 『TFT』 .
The mark of AC power in is 『AC 110V / 220V』 .

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
2.4 TO REPLACE FUSE


•Use screw driver to loose the Fuse corer.

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SECTION 3 : OPERATING INTRODUCTION

3.1 CONTROL PANEL

TriLite LED S600 series is equipped with light-head mounted control panel to control system on/off and intensity adjustment and to show intensity level status.

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
◦ TriLite LED S600 light-head is mounted with a film-keypad control panel for:


【Light-head Control Panel】



- (A) Power ON/OFF
- (B) Illumination Increase
- (C) Illumination Decrease
- (D) LED Failure Indicator
- (E) Illumination level Indicator
- (F) Endo. Light LED ON/OFF
- (G) ASM ON/OFF

70 Light-head Control Panel

50 Light-head Control Panel


- ⊙ Channel setting between remote control and lighthead has been completed before shipment. If user needs to reset, please follow the following procedure:
 1. Illumination of lighthead must be turned off.
 2. User needs to select a desired lighthead by pressing a button (L1、L2、or L3).
 3. Press and hold the power switch (A) of control box on lighthead.
 4. Aim the remote control to the signal receiver on lighthead and press the power switch (A) of remote control.
 5. The lighthead will be lit up to complete the resetting.
 6. ⊙ As soon as turn-on of lighthead power, the lighthead control will be turned-on. Users could select a desired lighthead by pressing a button (L1、L2、or L3) for illumination control. Press and hold power switch (A) for 1 second to turn-on the lighthead, while press and hold for 2 seconds to turn-off.
- ⊙ This remote control operates with power save mode, when the user stop operation for 10 seconds, CPU will enter into dormant mode until the user presses functional buttons(A、O、P、Q) to awaken CPU out of dormant mode and stay at the selected operation status prior to entering into dormant mode.
- ⊙ The remote control operation can be activated from a distance of up to 5 Meter radius.
- ⊙ The remote control uses infrared transmission to control the surgical light. Please avoid the presence of obstacles in the emission route of remote control for interfering.
- ⊙ The infrared transmission of the remote control is directional, please target the remote control to the direction of signal receiver Sensor.
- ⊙ Use two AAA batteries, make sure to take out batteries if the batteries are not used for a long time to avoid battery-liquid leakage and erosion on circuit board.
- ⊙ Battery life of this remote control is estimated to be three months of usage, make sure to replace battery periodically to ensure normal operation.
- ⊙ Please check the battery if the remote control does not work well.
- ⊙ The Camera control must be operated when any of O、P、Q is lit up.
- ⊙ When adjusting Endo. Light illumination, press F (Endo. Light), then select B and C to increase or decrease illumination.


3.2 Endo. Light

TriLite LED S600 has the fade lighting for MIS surgery, and to provide three LED lighting source in sterilizing handle, as shown in figure. The Endo. Light function be supported at Penetrating / Endo. Light control panel and remote controller. Push the function selector to star the Endo. Light LED and control the intensity of LED.

3.3 ASM (Automatic Shadow anagement only for TriLite LED S600-70)

TriLite LED S600-70 provides an option of ASM (Automatic Shadow Management) having 8 ASM sensors from 8 directions to detect any surgical team members' entrance into the illumination path range. ASM will automatically reduce the illumination for the obstructed LEDs and increase the illumination for the non-obstructed LEDs to deliver the excellent illumination on surgical site. ASM sensor is shown as the right picture.


3.4 INITIAL STATUS OF SYSTEM START

TriLite LED S600 series has the "automatic system self-turn-on design" whenever applied power is connected to the surgical light system to ensure turning the system itself on automatically as applied power coming back after temporary power shortage.

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


- By turning the sterilizable maneuvering handle to adjust size of the focused pattern or using optional focus control panel or remote control to adjust the focused pattern size.
- Users may reposition light-head through any positions with the sterilizable maneuvering handle or sterilizable handle.

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

PROBLEM	SOLUTION
Inner hub of suspension arm assembly drifts horizontally.	Shown as right figure, use Allen wrench M6 to tighten the setscrews, by turning clockwise. 
Outer hub (spring arm) of suspension arm assembly drifts horizontally.	Shown as right figure, use screwdriver to tighten the setscrews, by turning clockwise. 

Spring Arm drifts vertically.	1. Shown as right figure, insert and push Allen wrench M5 fully into the opening to adjust weight capacity. 2. To increase weight capacity, by turning Allen wrench toward counterclockwise. 3. To decrease weight capacity, by turning Allen wrench toward clockwise.	
Light-head movement height is too high or too low.	1. Shown as right figure, insert and push Allen wrench M5 fully into the opening to adjust spring arm height. Our factory set-up is 45°. 2. To have higher, fixed lowest position, by turning Allen wrench toward counter- clockwise. 3. To have lower, fixed lowest position, by turning Allen wrench toward clockwise.	

PROBLEM	SOLUTION	
Yoke contains rotation drift.	Shown as right figure, use bent nose plier to correct rotational drifting by turning the locknut inside the hub clockwise.	
Light-head contains rotational drift.	Shown as right figure, use bent nose plier to correct rotational drifting by turning the locknut inside the hub clockwise.	

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⚠ **WARNING:** Improper adjustment of spring arm could result in equipment damage.

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⚠ **WARNING:** Improper adjustment of spring arm could result in equipment damage.

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2. If problems still can't be solved, please feel free to contact the local distributor, or directly contact the Service Department to the **BenQ Medical**.

Mechanical		
Problem	Possible Cause	Suggested Correction
Sterilizable handle is not placed at correct position.	Incorrect setting of temperate and/or time during sterilization.	Inspect the positioning mechanism of sterilizable handle assembly.
	Structure becomes distorted because of ending of usage life.	Replace with a new sterilizable handle.
Lighthead Drafting	Ceiling structure unstable	Check ceiling structure
	Brake incorrectly adjusted	Adjust the brake
Lighthead too flexible or too rigid to manoeuvre	Insufficient lubrication	Call service department
	Brake incorrectly adjusted	Adjust the brake screw

Electricity		
Problem	Possible Cause	Suggested Correction
The lighthouse does not turn on	Improper installation	1. Use multi-meter to check if the voltage of main is correct for the input of power supply box. If no, check the main. 2. Check to have correct connection between anode and cathode from output of power supply box to the lamp. 3. Check if fuses in power supply box are burnt. If yes, replace with new fuses. 4. Use multi-meter to check power supply box to have 24 ± 0.5 Vdc. If no. Adjust DC voltage to reach 24 ± 0.5 Vdc. 5. Remove upper cover of lighthouse and check if the power cord is connected to the CN10 4-Pins connector on system board. If no. Connect them securely.
	Other reason	Call service department
Light flickers	Improper installation	Check the wiring in the lighthouse, power box.
	Component failure	Check conductive connector. in the carry arm
	Other reason	Call service department
Brightness decrease	Low voltage	Adjust the power supply unit output 24 ± 0.5 Vdc
	Component failure	1. Check conductive connector. in the carry arm. 2. Replace driver board. 3. Replace MC PCB 4. Replace driver board.
Nonresponded Control panel	Improper installation	Check the wiring between push button board and system board.
	Component failure	Replace push button board.
LED lifetime too short	Component failure	Replace LED MCPCB

Imaging		
Problem	Possible Cause	Suggested Correction
No image	No power supply for Camera	Check power supply
	No power supply for TFT-LCD	Check power supply of TFT-LCD.
	Incorrect or loose connection of video.	Check if connection of video wiring is correct.
	Parts failure.	1. Replace with a new driver board of camera. 2. Replace with a new camera
No focus	Object is not at center of screen.	Move camera to the center of screen.
	Status of manual focus adjustment	Manual focus adjustment (FAR · NEAR)
	No function for auto focus.	Switch AF/OFF to AF/ON.
	Parts failure	1. Replace with a new driver board of camera. 2. Replace with a new camera.
Blurred image	Lens is dirty.	Clean lens with lens-cleaning cloth
Image noise	Failure of conductive connector.	Call our technical department.

SECTION 5 : TECHNICAL DATA

5.1 TECHNICAL SPECIFICATIONS

LIGHTING DATA		TriLite LED S600-50	TriLite LED S600-70
Dia. of Light-head		56 cm (22.0")	75.6 cm (29.8")
Max. Intensity at 1m (39.4") Distance		Max. 120,000±10% Lux (11,148±10% Ft Candle)	Max. 160,000±10% Lux (14,864±10% Ft Candle)
Color Temperature		4,200±200 °K	42,00±200 °K
Color Rendering Index		95	95
Total Irradiance		370W/m ²	590W/m ²
Intensity Level Indication		12 Steps	12 Steps
Focused Pattern Adjustment	Minimum	16cm (6.3")	16 cm (6.3")
	Maximum	25cm (9.8")	31 cm (12.2")
Minimum Depth of Field without Refocusing		131 cm (51.6")	129 cm (50.8")
MECHANICAL DATA		TriLite LED S600-50	TriLite LED S600-70
Min. Height from Ground to Ceiling		300cm (118.1")	300cm (118.1")
Sterilizable Temperature of Maneuvering Handle		134 °C	134 °C
Light-head Weight		14.5 kg / 31.9 lb±5%	20.4 kg / 43.2 lb±5%
Total Weight	LED-50	48.0 kg±10% (105.7 lb)	
	LED-70	58.5 kg±10% (128.8 lb)	
	LED-75	92.8 kg ±10% (204.4 lb)	
ELECTRICAL DATA			
AC Input		100 - 240 V _{a.c.}	
Power Consumption		TriLite LED S600-50 : 191W	TriLite LED S600-70 : 288W
Number of LED		TriLite LED S600-50 : 96	TriLite LED S600-70 : 160
LED Color Rendering Index		95	
LED Average Life		≥ 30,000 Hours	
Applied Part		Type B • IEC 60601-1 • IEC 60601-2-41	
Fuses		5Ax2 、 10Ax1	
Batteries for Remote Control		AAA x 2	

◎ Built-in Endo. Lighting optical core assembly


ITEM	SPECIFICATION
Max. Intensity at 1m (39.4") Distance	300Lux ±20% Lux
Color Temperature	4,300±200 Kelvin
Power Consumption	15 W
Intensity Level Indication	12 steps/electron
LED Life	≥ 30,000 Hours

Operating/Shipping/Storage :


Operating/Shipping/Storage Environment:	
Ambient temperature	Shipping/Storage : -30 ~ 70°C
	Operating : 10 ~ 40°C
Relative humidity	Shipping/Storage : 20 ~ 80%
	Operating : 30 ~ 75%
Atmospheric pressure	Shipping/Storage : 500 ~ 1060 hPa
	Operating : 700 ~ 1060 hPa

5.2 RANGE OF CIRCULAR MOTION


© TriLite LED S600-755

© TriLite LED S600-750

© TriLite LED S600-700

TOTAL 29 PAGES (先 29 夏)

中華民國文件證明書
REPUBLIC OF CHINA TAIWAN
No. 1011000210249002
中華民國外交部
外交部文件證明書
This is to certify that the Signature/Seal of
Notary Public, Taipei District Court, is
authentic.
Date of Issue: Aug. 8, 2013
For the Minister and by authorization
Chu Tien-hsin
Secretary
Bureau of Consular Affairs
費用 fee: 公證費 NT\$400
G359790

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


TriLite LS800 Series Surgical Light

Installation/Operating Instructions

 BenQ Medical Technology Corporation
7F., No. 46, Zhou-Z St., Nei-Hu, Taipei 114, Taiwan

Customer Service e-mail :
Service@BenQMedicalTech.com
BenQMedicalTech.com
Customer Service Hotline :
Taiwan : 0800-015-533
China : +86 21-6327-7161-3 Ext. 812
International : +886 2-8797-5080 Ext. 5932

This operating instruction is intended to provide users with the proper information required for correct operation and maintenance of **BenQ** surgical light. In order to avoid danger to patient or users and/or damage to the equipment, please review the contents of this instruction thoroughly and carefully before use.

This information contained in this instruction is subject to change without prior notice.

BenQ STANDARD LIMITED WARRANTY

BenQ warrants all new surgical light sold by it directly or through a dealer or other authorized representatives, with expectation to replacement parts, spares, pads, and accessories, to be free from defects in material or workmanship, under normal use and service, for a period of one year. The warranty period shall begin with one year after shipment from **BenQ**.

BenQ's responsibility and liability shall be limited to the replacement of any parts which we, **BenQ**, determine to be defective within the applicable warranty period, labor will not be included in the warranty coverage. No dealer, agent, or employee is authorized to extend or enlarge this warranty.

BenQ Europe B.V.

Meerenakkerweg 1, 5652 AR, Eindhoven, Netherlands

TABLE OF CONTENTS

Summary of Safety and General Precautions	_____ Page 2
Section 1 : Section Introduction	
1.1 General	_____ Page 5
1.2 Contents of Wooden Case	_____ Page 6
1.3 Power Source Requirement	_____ Page 6
1.4 Explanation of Precaution Labels	_____ Page 7
1.5 Cleaning or Disinfecting	_____ Page 10
Section 2 : Installation Introduction	
2.1 Ceiling Mounting Methods	_____ Page 12
2.2 Installation Procedure of Suspension Arm and Light-head	_____ Page 15
2.3 Power Connection Instructions for Ceiling Power Box	_____ Page 18
2.4 To replace fuse	_____ Page 19
2.5 The specification of fuse	_____ Page 19
Section 3 : Operating Introduction	
3.1 Control Panel	_____ Page 20
3.2 Endo. Light Control	_____ Page 22
3.3 Initial Status of System Start	_____ Page 22
3.4 Focused Pattern Size Adjustment and Light-head Positioning	_____ Page 22
3.5 Brake/Friction Adjustment	_____ Page 23
3.6 Spring Arm Adjustment	_____ Page 25
Section 4 : Adjustment and Maintenance	
4.1 Troubleshooting	_____ Page 26
Section 5 : Technical Data	
5.1 Technical Specifications	_____ Page 29
5.2 Range of Circular Motion	_____ Page 33
5.3 Operating/Shipping/Storage Environment	_____ Page 34
Safety Declaration	_____ Page 35

SUMMARY OF SAFETY PRECAUTIONS

The user has the ultimate responsibility of preventing damage to the light and nearby equipment or injury to the patient and surgical team.

Warning: indicating possibility of personal injury.

Caution: indicating potential damage to equipment.

Note: indicating pertinent facts or conditions.

WARNING:

PINCHING HAZARD:

- ⊙ Pinch points are composed in extreme light-head articulations.
- ⊙ User should be aware of and examine these points before operating the light.
- ⊙ Potential pinch and/or crush hazard exists if hands or fingers are placed between edges of yoke and light-head during extreme spring arm posturing.

COLLISION HAZARD:

- ⊙ Please put the light to an appropriate place after use to avoid collision.

FAILURE HAZARD:

- ⊙ Do not use **BenQ** light parts for other than their stated function or parts manufactured by other manufacturers to avoid possible patient injury, user injury, or surgical light or accessory failure. Any warranties could also be rendered void by the use of non-original parts.
- ⊙ For safe and reliable operation of the lamp, installation, maintenance and repair of the equipment must be done only by qualified service personnel. Preventive maintenance of the lamp requires a regular schedule.
- ⊙ Serious injury or lamp damage may result if user of the lamp is not fully familiar with the controls for lamp operation and needed positioning.

ELECTRICAL SHOCK HAZARD:

- ⊙ Improper electrical grounding or no grounding may result in personal injury or equipment damage.
- ⊙ Do not remove cover of electrical box and/or light-head before disconnecting electrical power to the system or other than those stated in this instruction.
- ⊙ Do not attempt to replace bulb unless power is turned off and light-head is cool adequately.
- ⊙ Do not attempt to clean light-head unless power is turned off and light-head is cool adequately.

EXPLOSION HAZARD:

Do not use in the presence of flammable anesthetics.

BURNT HAZARD:

- ⊙ The total irradiance of some series of this product might be more than 1,000 W/m², please be attention to this effect.

ELECTROMAGNETIC INTERFERENCE HAZARD:

- ⊙ The light has been approved by IEC60601-1 and IEC60601-1-2 of Medical electrical equipment standard.
- ⊙ Do not use the light near the equipment that may cause electromagnetic interference, e.g. mobile phone, radio reception equipment, MRI.
- ⊙ Observe the anti-electromagnetic interference of EMC

HIGH-FREQUENCY EQUIPMENT CONTACT HAZARD:

- ⊙ Contact with conductivity when using high-frequency equipment such as high frequency electrocautery, defibrillator, defibrillator monitor may cause personal burning. Personal should not contact with metal component or wet surgical cloth of high-frequency equipment.

WASTE HAZARD:

- ⊙ Certain component should be disposed by the company with hazard waste disposal licence, and should comply with local regulations

PERSONAL INJURY HAZARD:

- ⊙ Do not remove or modify light fixture by unauthorized technicians. Refer servicing to authorized technicians.
- ⊙ Do not operate the light if the Heatsorb filters or cover lens are broken or removed. Operation without the Heatsorb filters or cover lens or with broken Heatsorb filters or cover lens could result in injury to the patient or user.

CAUTION:

EQUIPMENT DAMAGE:

1. Plug into proper socket and confirm grounding.
2. The surgical light doesn't install a disconnecter. Power supply is in the normal open state. Connect each power supply input to a disconnecter, or a switch with circuit break or knife switch.
3. Prevent leakage of fluids into inside of light-head.
4. Do not expose the lamp to excessive moisture.
5. Do not use harsh cleaning, disinfecting, or anti-static liquid on the light. Use of phenolic-, iodophor-, glutaraldehyde-, alcohol-based disinfectants on surface of the light could cause some degree of staining, pitting, and/or discoloration damage on painting or optical cover.
6. Use only clean, soft, and lint-free cloth to wipe both internal or external surfaces of this light when cleaning.
7. Adjust to the lowest intensity before turn down the light.
8. Do not illuminate front lens of surgical light by UV sterilization light in order to prevent the aging of front lens.
9. Use original component.

WARNING:

1. Electrical devices (e.g. cell phones, radio, magnetic resonance tomographs) can interfere with the functioning of the product when used near the product. Please note the information concerning electromagnetic compatibility (EMC) (radiation and resistance to interference) contained in the Safety Declaration. Adhere to those specifications when using electrical devices and respond properly in the event of effects on the device or the product.
2. Switch off your wireless device when requested to do so in hospitals, clinics or health care facilities. These requests are designed to prevent possible interference with sensitive medical equipment.

Listed below are all the accessories that may be attached to the surgical light.
Power cord / 750mm $\pm$ 50 mm

Figure Illustration :



	Warning: indicating possibility of personal injury.
Model	Model
REF	Product Number
SN	Serial Number
	Grounding
	Neutral connection point of permanently installed equipment
	Power On/Off
	S.M.S. On/Off
	Illumination Decrease
	Illumination Increase
	Endo. Light LED On/Off
	CCD zoom in
	CCD zoom out
	CCD CW Rotate
	CCD CCW Rotate

SECTION 1 : SECTION INTRODUCTION

1.1 GENERAL

BenQ's TriLite LS800 series surgical lights use LED technology to provide a field of clear illumination in surgical, diagnosis, or treatment applications. The single light-head model is intended for general and large surgical, diagnosis, and/or treatment applications and may be used in the situation, where the optimal illumination, cold light, color temperature and color rendering index suitable for procedures observations can be provided where blood vessels, organ tissues, and nerves of surgical site can be differentiated. The color temperature of TriLite LS800 is $4,300\pm 200^{\circ}\text{K}$ that is close to average day light, operating rooms for all surgical procedures and where a safe condition is required.

This lighting condition will come in handy to all surgeons as if all the surgical procedures are performed at 8 o'clock in the morning of summertime. TriLite LS800 surgical lights series deliver cold light required for surgical procedures where blood coagulation and tissue dryness must be prevented and provide surgical team with good working environment and eliminate chances of uncomfortable illumination owing to overheat.

1.2 CONTENTS OF CARTON CASE

There are cartons for suspension arm assembly and light head assembly separately inside the carton case.

Accessories:

1. Suspension Arm Assembly - not more than 3 pcs.;
2. Spring Arm - not more than 3 pcs.;
3. Mounting plate - 1 pcs.;
4. Ceiling housing - 1 pcs.;
5. Mounting kit - 1 pcs.;
6. Power Box - 1 pcs.;
7. Light-heads - 1 pcs.;
8. Lighting unit with built-in CCD - 1 pcs.;
9. Optical Core Handle Assembly (with Sterilizable Handle) - 1 pcs.;
10. Remote control - 1 pcs.;
11. Separate CCD - 1 pcs.;
12. Bracket for mounting a camera - not more than 1 pcs.;
13. TFT-display - not more than 3 pcs.;
14. Control unit wall - 1 pcs.;
15. Medical PC ORLIS - 1 pcs.;
16. Handle with built-in LEDs EndoLED - 1 pcs.;
17. Vertical Gimbals Joint - not more than 3 pcs.;
18. LEDs - not more than 160 pcs.;
19. AAA batteries - not more than 2 pcs.;
20. Power cable - not more than 5 pcs.;
21. Instruction for use - 1 pcs.;
22. Extension pipe - not more than 2 pcs.

1.3 POWER SOURCE REQUIREMENT

1. TriLite LS800 series surgical light can accept applied power source 110/220Vdc, 50/60Hz. (Please see the label on the power box before connecting power source.)
2. If there is a need to extend the length of power cable between the power box and surgical light, please follow the following formula to calculate the precise cross section area for each individual core.

$$A=L(m) \times P(w) \times 1/900 \quad (\text{mm}^2/\text{Wm})$$

Note: The cross section area of the chosen power cable must greater than "A" value.

A= The cross section area of the extended power cable (mm²). L= The length of extended power cable between power box and surgical light (m). P= Power consumption (W), please refer the power consumption to the light-head label.

⚠CAUTION:

Always limit the applied power source at the range of ±10% to avoid equipment damage.

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

【LABEL ATTACHED PLACE AND EXPLANATION】

Note: All of "X" mark on each label is only for explanation purpose. It will be varied as serial number for each different unit.

1.5 CLEANING/DISINFECTION/STERILIZATION

Users must contact hospital's sanitary specialists. The recommended products and procedures must be applied. Should there be any doubt concerning the compatibility of active agents to be used contact our local distributor or the local BenQ Medical Customer Service Center.

A. Cleaning and disinfecting the surgical light

1. Cleaning and disinfecting :

- * Check that the power is switched off.
- * Remove the sterilizable handles.
- * Clean the system using a cloth soaked in surface detergent, in line with the manufacture's recommended dilution, application time and temperature.
- * Use a cloth to rinse the unit with clean water and wipe dry.
- * Disinfect the system using a cloth soaked in disinfectant, in line with the manufacture's recommendations.
- * Use a cloth to rinse with clean water in order to removes residues (in particular products containing aldehydes, quaternary ammonium or surfactants).
- * Wipe with a dry cloth .

2. Recommended products:

Getinge USA: TEC-QUAT 256.

Anios: SURFA'SAFE , ANIOSYME DLT PLUS

3. Examples of prohibited products

Solutions containing glutaraldehyde, phenol, iodine, bleach, alcohol or chloride ions must not be used.

Fumigation methods are unsuitable for disinfecting the unit and must not be used.

B. Cleaning and sterilizing the handles

1. Before cleaning

- * Use a soft cloth immediately after use to wipe away soiling from the handle surface.
- * Store handles in a place that keeps them moist to make further cleaning easier.

2. Cleaning

- * Soak the handles in a detergent solution. (note.1)
- * Soak for 15 minutes to allow the solution to act, then clean by hand with a soft brush and a lint-free cloth.
- * During cleaning, check regularly that the handles are fully clean and that no soiling remains on the inside or outside.
- * If any soiling remains, repeat cleaning or use an ultrasonic cleaning process.
- * Rinsing: Rinse thoroughly in clean water to completely eliminate the detergent solution.
- * Drying: Wipe with a clean lint-free cloth.

3. Disinfection and Sterilization

Handles(note. 2) may be disinfected by clean machine(Ultrasonic) , in line with the manufacture's recommendations.

After cleaning, the handles must be steam sterilized as set out below:

Countries	Sterilization Cycle	Temperature [°C]	Time[min]	Drying[min]
Taiwan	Prevacuum	132~135	10	16
USA	Prevacuum	132~135	10	16
Other	Prevacuum	Comply with national regulations		

- * Check that each handle is clean before continuing the process.
- * Wrap the handles with sterilization wrapper material (double wrapper or equivalent).
Handles may also be placed in paper or plastic sterilization bags (note.3), for easier identification and reuse.
- * Place the handles on sterilizer trays with the opening downwards.
- * Package with biological and/or chemical indicators for monitoring the sterilization process, in accordance with applicable regulations.
- * Run the sterilization cycle according to the sterilizer manufacturer's instructions.
 - To ensure correct sterilization do not allow any soiling to penetrate inside the handle.
 - Handles are not guaranteed beyond 350 sterilization cycles with the above sterilization parameters.
- * Dispose of sterilizable handles in the same way as other hazardous products in a hospital environment.

note.1: Never soak the handles in enzyme-based detergents as they may damage the material used; rinse thoroughly if these detergents are used.

note.2: The material of handles is Polysulfone (PSF).

note.3: Sterilization bag suppliers: Sigma Medical Supplies Corp · Medical Action Industries.

2.1 CEILING MOUNTING METHODS

Existing Ceiling Structure

Finished Ceiling

Finished Floor

Figure (2.1.1)

Existing Ceiling Structure

Mounting Disc Plate


Nuts

10

47

Figure (2.1.2)

Mounting Disc Plate

13

- ✓ If a ceiling height exceeds 9'10"(300cm), a properly reinforced intermedium frame is required to support the lamp. If there is a finished ceiling, a properly reinforced intermedium frame is also needed for installation of the lamp. The following is the methods of a properly reinforced intermedium frame, depending on ceiling height.

Situation (2) If the distance between existing ceiling structure and finished ceiling is within 11.8"(30cm), the installation instruction of a properly reinforced intermedium frame shown as figure (2.1.3).

Use expansion (foundation) bolts to fix the 4 holes, $\varnothing 0.87"$ (22mm) of mounting disc plate in place, then use 4 connect nuts and studs to extend your needed length. Ceiling plate must be level in two directions.

The upper surface of the mounting disc plate must align with the finished floor shown as figure (2.1.3).


Figure (2.1.3)

- 1: Put M16 expansion (foundation) bolts into existing ceiling structure for 3.15"(8cm).
- 2: 3/4"-10UNC studs(or M18) studs.

Situation (3) If the distance between existing ceiling structure and finished ceiling is from 11.8" (30cm) to 23.6" (60cm), refer to figure (2.1.4).

Base on installation method in the situation (2), add a reinforced tube on outside of each of the 4 expansion (foundation) bolts and studs to reinforce the whole structure.


Figure (2.1.4)

Reinforced tubes:
 $\varnothing 1.22"$ (31mm, inner) $\times \varnothing 1.5"$ (38mm, outer).

Situation (4) If the distance between existing ceiling structure and finished ceiling is from 23.6" (60cm) to 47.24" (120cm), refer to figure (2.1.5).

If distance between existing ceiling structure and finished ceiling is 39.4" (100cm), prepare a 80cm of height of large-tube extension frame and save 20cm of height for incorporating with the situation (1) method to connect large tube with mounting disc plate.

This will effectively adjust length and level the whole reinforced intermedium frame.


Figure (2.1.5)

Warning: It must be connected with grounding conductor to protect the circuit of surgical light.

✓ **Manufacturing of extension tube frame:**

- Length of extension tube between each point on both cutting edges must be the same. This will ensure that frame is always align and level after weld with the steel plates on each end shown as figure (2.1.6).









Figure (2.1.6)


Warning: It must be connected with grounding conductor when you install mounting disc plate and extension tube.

2.2 INSTALLATION PROCEDURE OF SUSPENSION ARM AND LIGHT-HEAD


1  Mounting Disc Plate Locking Bolt Suspension Arm Assembly • Fix suspension arm assembly onto mounting disc plate with 6 locking bolts. Note: Always level mounting disc plate before fixing suspension arm assembly.	2  • Remove the Set Screw on the Mount of Spring Arm.
3  • Insert the Male Connector into Mount on Spring Arm, and fixed the Male Connector with Set Screw.	4  • Insert the Mount of Spring Arm onto Extension Arm.
5  • Align braking setscrew hole with a concavity on inner plastic tube precisely and screw a braking screw to have proper brake on the spring arm. Caution: Always align braking setscrew hole with a concavity on inner plastic tube precisely before screwing a braking setscrew.	6  • Use screwdriver to properly tighten braking setscrew, and do micro-adjustment after light-head is installed.


- Remove screw on spring arm with phillips screw driver.


- Lift safety sleeve slightly, twist to the right, and then pull right up.
- Remove safety segment.

- Insert spigot of light-head yoke into socket on spring arm.




- Replace safety segment back to socket.
- Remove protection sheet of cover lens of light-head.

- Pull the safety sleeve downwards.
- Screw the locking screw back to lock the sleeve with the socket.

- Attaching this label to spring arm sleeve.
- Note :** Only qualified service personnel is allowed to peel off the label once label has been applied.

13  Ceiling Shroud • Twist slightly and place ceiling shroud to cover central base of suspension arm assembly.	14  • Twist and turn the ceiling shroud again to pass upper suspension arm.
15  Setscrew • Lock ceiling shroud with setscrews.	

2.3 POWER CONNECTION INSTRUCTIONS FOR CEILING POWER BOX


Figure (2.3.1) Power Box (Switching Power Supply)

(I)
Suspension Arm Identification
Label-Lamp(1)-Lamp(1)-Red

(Red)

(J)
Suspension Arm Identification
Label-Lamp(2)-Lamp(2)-Green

(Green)

(K)
Suspension Arm Identification
Label-Lamp(3)-Lamp(3)-Blue

(Blue)


The mark of primary light head is 『Lamp (1)』 .
The mark of secondary light head is 『Lamp (2)』 .
The mark of third light head is 『Lamp (3)』 .
The mark of the separate carrying arm of TFT is 『TFT』 .
The mark of AC power in is 『100-240Va.c.』 .

⚠ WARNING:

- Improper electrical grounding or no grounding may result in personal injury or equipment damage.
- To check the mark, before connecting with V_{DC} and power cable.
- Each light head has its own power supply box. This means that dual light head surgical light has 2 power supply boxes.


2.4 TO REPLACE FUSE


•Use screw driver to loose the Fuse cover.

•Replace the burn out Fuse with the same specific Fuse.


•Reinstall the Fuse cover.

•Use screw driver to lock the Fuse cover and then finish the Procedure.

2.5 THE SPECIFICATION OF POWER SUPPLY AND FUSE

Fuse specification :
Fuse1&Fuse2= 5A/250V Slow blow, Glass tube type.

SECTION 3 : OPERATING INTRODUCTION


3.1 CONTROL PANEL


TriLite LS800 series is equipped with light-head mounted control panel to control system on/off and intensity adjustment and to show intensity level status.

※The function buttons on the light-head control panel have internationally recognized symbols and function name for all functions to allow easy identification.

【Light-head Control Panel】

- (A) Power ON/OFF
- (B) Brightness Increase
- (C) Brightness Decrease
- (D) Brightness Level Indicator
- (E) Endo. Light LED ON/OFF
- (F) S.M.S. ON/OFF
- (G) CCD zoom in
- (H) CCD zoom out
- (I) CCD CW rotate
- (J) CCD CCW rotate


- ⊙ Channel setting between remote control and lighthouse has been completed before shipment. If user needs to reset, please follow the following procedure:
 1. Illumination of lighthouse must be turned off.
 2. User needs to select a desired lighthouse by pressing a button (L1、L2、or L3).
 3. Press and hold the power switch (A) of control box on lighthouse.
 4. Aim the remote control to the signal receiver on lighthouse and press the power switch (A) of remote control.
 5. The lighthouse will be lit up to complete the resetting.
 6. ⊙ As soon as turn-on of lighthouse power, the lighthouse control will be turned-on. Users could select a desired lighthouse by pressing a button (L1、L2、or L3) for illumination control. Press and hold power switch (A) for 1 second to turn-on the lighthouse, while press and hold for 2 seconds to turn-off.
- ⊙ This remote control operates with power save mode, when the user stop operation for 10 seconds, CPU will enter into dormant mode until the user presses functional buttons(A、O、P、Q) to awaken CPU out of dormant mode and stay at the selected operation status prior to entering into dormant mode.
- ⊙ The remote control operation can be activated from a distance of up to 5 Meter radius.
- ⊙ The remote control uses infrared transmission to control the surgical light. Please avoid the presence of obstacles in the emission route of remote control for interfering.
- ⊙ The infrared transmission of the remote control is directional, please target the remote control to the direction of signal receiver Sensor.
- ⊙ Use two AAA batteries, make sure to take out batteries if the batteries are not used for a long time to avoid battery-liquid leakage and erosion on circuit board.
- ⊙ Battery life of this remote control is estimated to be three months of usage, make sure to replace battery periodically to ensure normal operation.
- ⊙ Please check the battery if the remote control does not work well.
- ⊙ The Camera control must be operated when any of O、P、Q is lit up.
- ⊙ When adjusting Endo. Light illumination, press F (Endo. Light), then select B and C to increase or decrease illumination.

3.2 Endo. Light

TriLite LS800 has the fade lighting for MIS surgery, and to provide 15 LED lighting source around sterilizing handle, as shown in figure.

The Endo. Light function is supported at Endo. Light control panel and remote controller.

Push the function selector to start the Endo. Light LED and control the intensity of LED.


3.3 INITIAL STATUS OF SYSTEM START

TriLite LS800 series has the "automatic system self-turn-on design" whenever applied power is connected to the surgical light system to ensure turning the system itself on automatically as applied power coming back after temporary power shortage.


3.4 FOCUSED PATTERN SIZE ADJUSTMENT AND LIGHT-HEAD POSITIONING



- ⊙ By turning the sterilizable maneuvering handle to adjust size of the focused pattern or using optional focus control panel or remote control to adjust the focused pattern size.
- ⊙ Users may reposition light-head through any positions with the sterilizable maneuvering handle or sterilizable handle.

3.5 BRAKE/FRICTION ADJUSTMENT

- Use operating instructions to verify trouble symptoms.
- After symptom has been verified, refer to the following table and select example that is most appropriate to your problem. Follow recommended correction.

PROBLEM	SOLUTION
Inner hub of suspension arm assembly drifts horizontally.	Shown as right figure, use Allen wrench M6 to tighten the setscrews, by turning clockwise. 
Outer hub (spring arm) of suspension arm assembly drifts horizontally.	Shown as right figure, use screwdriver to tighten the setscrews, by turning clockwise. 
Spring Arm drifts vertically.	1. Shown as right figure, insert and push Allen wrench M5 fully into the opening to adjust weight capacity. 2. To increase weight capacity, by turning Allen wrench toward counterclockwise. 3. To decrease weight capacity, by turning Allen wrench toward clockwise. 
Light-head movement height is too high or too low.	1. Shown as right figure, insert and push Allen wrench M5 fully into the opening to adjust spring arm height. Our factory set-up is 45°. 2. To have higher, fixed lowest position, by turning Allen wrench toward counter- clockwise. 3. To have lower, fixed lowest position, by turning Allen wrench toward clockwise. 

Yoke contains rotation drift.	Shown as right figure, use bent nose plier to correct rotational drifting by turning the locknut inside the hub clockwise. 
Light-head contains rotational drift.	Shown as right figure, use bent nose plier to correct rotational drifting by turning the locknut inside the hub clockwise. 

※If you still can not solve problem, please contact our local agent or our Customer Service Center in Taiwan with detailed problem, condition, and series number of a problem unit.

WARNING: Improper adjustment of spring arm could result in equipment damage.

3.6 SPRING ARM ADJUSTMENT

A. Height movement adjustment:

The downward movement of spring arm is restricted at a fixed lowest position. If you need a further downward position beyond the fixed lowest position, you may adjust with an Allen wrench.

1. Insert and push Allen wrench M5 fully into the opening **X** shown as figure(a).
2. To have **higher** fixed lowest position, by turning Allen wrench toward **counter-clockwise**.
3. To have **lower** fixed lowest position, by turning Allen wrench toward **clockwise**.

B. Weight bearing capacity adjustment:

There is only 1 point of adjustment to correct any downward drifting or weight bearing capacity of the spring arm.

1. Insert and push Allen wrench M5 fully into the opening **Y** shown as figure(a).
2. To **increase** weight bearing capacity, by turning Allen wrench toward **counterclockwise**.
3. To **decrease** weight bearing capacity, by turning Allen wrench toward **clockwise**.

 **WARNING:** Improper adjustment of spring arm could result in equipment damage.


Figure (a) Spring Arm Adjustment

SECTION 4 : ADJUSTMENT AND MAINTENANCE

4.1 TROUBLESHOOTING

1. Use operating instructions to verify trouble issue, after symptom has been verified, refer to Exploded Views and Parts Lists and the following table and select example that is most appropriate to your problem. Follow recommended correction for guidance when attempting repairs.

2. If problems still can't be solved, please feel free to contact the local distributor, or directly contact the Customer Service Center to the **BenQ Medical**.

Mechanical		
Problem	Possible Cause	Suggested Correction
Sterilizable handle is not placed at correct position.	Incorrect setting of temperate and/or time during sterilization.	Inspect the positioning mechanism of sterilizable handle assembly.
	Structure becomes distorted because of ending of usage life.	Replace with a new sterilizable handle.
Lighthead Drafting	Ceiling structure unstable.	Check ceiling structure.
	Brake incorrectly adjusted.	Adjust the brake.
Lighthead too flexible or too rigid to manoeuvre	Insufficient lubrication.	Please call our local distributor or customer service center.
	Brake incorrectly adjusted.	Adjust the brake screw.

24

Electricity		
Problem	Possible Cause	Suggested Correction
The lighthouse does not turn on	Improper installation.	1. Use multi-meter to check if the voltage of main is correct for the input of power supply box. If no, check the main. 2. Check to have correct connection between anode and cathode from output of power supply box to the lamp. 3. Check if fuses in power supply box are burnt. If yes, replace with new fuses. 4. Use multi-meter to check power supply box to have $24 \pm 0.5 V_{d.c.}$. If no. Adjust DC voltage to reach $24 \pm 0.5 V_{d.c.}$. 5. Remove upper cover of lighthouse and check if the power cord is connected to the CN10 4-Pins connector on system board. If no. Connect them securely.
	Other reason.	Please call our local distributor or customer service center.
Light flickers	Improper installation.	Check the wiring in the lighthouse, power box.
	Component failure.	Check conductive connector. in the carry arm.
	Other reason.	Please call our local distributor or customer service center.
Brightness decrease	Low voltage	Adjust the power supply unit output $24 \pm 0.5 V_{d.c.}$
	Component failure.	1. Check conductive connector. in the carry arm. 2. Replace driver board. 3. Replace MC PCB. 4. Replace driver board.
Nonresponded Control panel	Improper installation.	Check the wiring between push button board and system board.
	Component failure.	Replace push button board.
LED lifetime too short	Component failure.	Replace LED MCPCB.

Imaging		
Problem	Possible Cause	Suggested Correction
No image	No power supply for Camera.	Check power supply.
	No power supply for TFT-LCD.	Check power supply of TFT-LCD.
	Incorrect or loose connection of video.	Check if connection of video wiring is correct.
	Parts failure.	1. Replace with a new driver board of camera. 2. Replace with a new camera.
No focus	Object is not at center of screen.	Move camera to the center of screen.
	Status of manual focus adjustment.	Manual focus adjustment (FAR , NEAR).
	No function for auto focus.	Switch AF/OFF to AF/ON.
	Parts failure.	1. Replace with a new driver board of camera. 2. Replace with a new camera.
Blurred image	Lens is dirty.	Clean lens with lens-cleaning cloth.
Image noise	Failure of conductive connector.	Please call our local distributor or customer service center.

SECTION 5 : TECHNICAL DATA

5.1 TECHNICAL SPECIFICATIONS

LIGHTING DATA		TriLite LS800-500	TriLite LS800-700
Dia. of Light-head		64 cm (25.2")	76.4 cm (30")
Max. Intensity at 1m (39.4") Distance		120,000 Lux	160,000 Lux
Color Temperature		4,300±200°K	4,300±200°K
Color Rendering Index		≥95	≥95
Total Irradiance		305 W/m ²	421 W/m ²
Intensity Level Indication		12 Steps	12 Steps
Light field diameter	Light Field d50	16 cm	18 cm
	Light Field d10	28 cm	32 cm
Depth of illumination (L1+L2)		130 cm	120 cm
	Remaining illuminance with one mask	45%	79%
	Remaining illuminance with two masks	47%	52%
	Remaining illuminance with tube	95%	85%
	Remaining illuminance with tube and one mask	45%	67%
	Remaining illuminance with tube and two masks	44%	47%
MECHANICAL DATA		TriLite LS800-500	TriLite LS800-700
Min. Height from Ground to Ceiling		290 cm	290 cm
Light-head Weight		14.5 kg / 31.9 lb±5%	18.5 kg / 40.7 lb±5%
ELECTRICAL DATA			
AC Input		100 - 240 V _{a.c.}	
Power Consumption		TriLite LS800-500 : 95W	TriLite LS800-700 : 135W
Number of LED		TriLite LS800-500 : 57	TriLite LS800-700 : 95
LED Average Life		≥40,000 Hours(in the best situation)	
Applied Part		IEC 60601-1、IEC 60601-2-41	
Fuses		5Ax2、10Ax1	
Battery of Remote Control		2 AAA batteries	

※Test Temperature Standard (20 °C)

30


(HD CCD Specifications)

Items		Specifications
Style No.		Sony FCB-H11
Signal system	HD	1080/59.94i, 1080/50i, 720/59.94p, 720/50p
	SD	NTSC (Crop/Squeeze), PAL (Crop/Squeeze)
Image sensor		1/3-type CMOS
Resolution		2,000,000 pixels
Minimum illumination		1 lx
Video output	HD	Analog component : Y/Pb/Pr
	SD	VBS : 1.0 Vp-p (sync negative) Y/C
White balance		2500-8000 k
Lens	Optical	×10 F=5.1mm(wide) to 51.0mm(tele)
	Digital	×12 (120X with optical zoom)
Power requirements		6-12 Vd.c.
Power consumption		4.8W
Dimensions		47.2(W)×43.1(H)×72.2(D) mm
Mass		120 g

(TFT LCD Monitor)


Items		Specifications
Style No.		SONY LMD-2110
Resolution		1920×1080 pixels
Colors		Approx. 16.7 million colors
Brightness		600 cd/m ² (DICOM preset at 400 cd/m ²)
Input signal / terminal	Composite	BNC(×1), 1.0Vp ± 3db sync negative
	Y/C	Mini DIN 4-pin(×1) Y: 1.0Vp ± 3db sync negative C: 0.286Vp-p ± 3db
	RGB	BNC×3 RGB: 0.7Vp-p ± 3db
	HDMI	HDMI×1
Viewing angle		170°(H) × 160°(V)
Power requirements		AC 100V to 240V, 50/60Hz
Operating temperature		0°C to 35°C ; Recommended: 20°C to 30°C
Operating humidity		30% to 85% (no condensation)
Power consumption		Maximum approx. 56 W, 1.0 A to 0.5 A
Dimensions		515(L)×86 (W)×355(H) mm
Mass		6.9 kg

// Front Panel

- ① Tally lamp ② standby switch and Indicator
- ③ KEY INHIBIT Indicator ④ VOLUME buttons
- ⑤ Menu operation buttons ⑥ ASPECT select button
- ⑦ SCAN select button ⑧ BLUE ONLY button
- ⑨ EXT SYNC (external sync) button ⑩ SDI button
- ⑪ RGB/COMPONENT button ⑫ HDMI button ⑬ LINE

// Connector Panel

- ① HDMI IN connector ② HDMI cable holder
- ③ OPTION In connector ④ OPTION AUDIO In (Phono Jack)
- ⑤ EXT SYNC In/Out (external sync) (BNC)
- ⑥ RGB/COMPONENT (BNC), Audio (Phono Jack)
- ⑦ PARALLEL REMOTE (modular connector)
- ⑧ LINE [composite (BNC), Y/C (Mini DIN 4-pin), Audio (Phono Jack)]
- ⑨ AC In


(TFT Mounting Bracket)

Items	Specifications
The maximum rotation angle of the mounting bracket on the spring arm axis	$\geq 280^\circ$
maximum action diameter of TFT mounting bracket	338.8cm
Maximum pitch axis tilt of TFT mounting bracket	$195^\circ \pm 5^\circ$
Voltage	100-240 Va.c. · 50/60Hz


Model Series	
Model	TriLite LS800
Series Style No.	TriLite LS800-500X1X2X3X4 · TriLite LS800-700 X1X2X3X4 · TriLite LS800-550 X1X2X3X4 · TriLite LS800-770 X1X2X3X4 · TriLite LS800-750 X1X2X3X4 · TriLite LS800-555 X1X2X3X4 · TriLite LS800-755 X1X2X3X4 X1X2X3X4 is optional : H (Built-in CCD Camera) C (External CCD Camera) T (TFT display system) TT (Dual TFT display system) W (Wall control box system) S (Snake fiber optic assembly) M (Mobile Stand system) (For TriLite LS800-500 only) B (Battery system)

5.2 Range of Circular Motion

© TriLite LS800-75H (75,55,7HT)

•TriLite LS800-75HT (755H,77HT,775H, 755H,755,555)

5.3 Operating/Shipping/Storage Environment

Operating/Shipping/Storage Environment:	
Ambient temperature	Shipping/Storage : -20 ~ 50°C
	Operating : 10 ~ 40°C
Relative humidity	Shipping/Storage : 20 ~ 80%
	Operating : 30 ~ 75%
Atmospheric pressure	Shipping/Storage : 500 ~ 1060 hPa
	Operating : 700 ~ 1060 hPa

Safety Declaration

Electromagnetic emissions and Immunity for TriLite LS800 series (110V a.c. models)

Recommended separation distances between portable and mobile RF communications equipment and the ME equipment			
The TriLite LS800 (AC 110V 60Hz) is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the TriLite LS800 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the TriLite LS800 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{F_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.7	3.7	7.37
100	11.67	11.67	23.33

**Declaration – electromagnetic emissions and immunity –
for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING and
are specified for use only in a shielded location**

are specified for use only in a shielded location

The TriLite LS800 declaration – electromagnetic immunity			
The TriLite LS800 system is intended for use in the electromagnetic environment specified below. The customer or the user of the TriLite LS800 system should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3V	Portable and mobile RF communications equipment should be used no closer to any part of the EQUIPMENT or SYSTEM including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Interference may occur in the vicinity of equipment marked with the following symbol. 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m	

Declaration – electromagnetic immunity			
The TriLite LS800 system is intended for use in the electromagnetic environment specified below. The customer or the user of the TriLite LS800 system should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the EQUIPMENT or SYSTEM requires continued operation during power mains interruptions, it is recommended that the EQUIPMENT or SYSTEM be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Declaration – electromagnetic emissions		
The TriLite LS800 is intended for use in the electromagnetic environment specified below. The customer or the user of the TriLite LS800 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
CE emissions CISPR11	Group 1	The TriLite LS800 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The TriLite LS800 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RE emissions CISPR11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Complies	

Electromagnetic emissions and Immunity for TriLite LS800 series (230Va.c. models)

Recommended separation distances between portable and mobile RF communications equipment and the ME equipment			
The TriLite LS800 (AC 230V 50Hz) is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the TriLite LS800 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the TriLite LS800 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{f_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.7	3.7	7.37
100	11.67	11.67	23.33

**Declaration – electromagnetic emissions and immunity –
for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING and
are specified for use only in a shielded location**

The TriLite LS800 declaration – electromagnetic immunity			
The TriLite LS800 system is intended for use in the electromagnetic environment specified below. The customer or the user of the TriLite LS800 system should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3V	Portable and mobile RF communications equipment should be used no closer to any part of the EQUIPMENT or SYSTEM including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Interference may occur in the vicinity of equipment marked with the following symbol. ((())) 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m	

Declaration – electromagnetic immunity			
The TriLite LS800 system is intended for use in the electromagnetic environment specified below. The customer or the user of the TriLite LS800 system should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the EQUIPMENT or SYSTEM requires continued operation during power mains interruptions, it is recommended that the EQUIPMENT or SYSTEM be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Declaration – electromagnetic emissions		
The TriLite LS800 is intended for use in the electromagnetic environment specified below. The customer or the user of the TriLite LS800 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
CE emissions CISPR11	Group 1	The TriLite LS800 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The TriLite LS800 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RE emissions CISPR11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Complies	

TOTAL 30 PAGES (先 39 夏)

中華民國文件證明書
REPUBLIC OF CHINA (TAIWAN)
No. 1041000023041
中華民國 102 年 08 月 08 日
MINISTRY OF FOREIGN AFFAIRS
REPUBLIC OF CHINA
This is to certify that the Signature/Seal of
Notary Public, Taipei District Court, is
authentic.
Date of issue: Aug. 8, 2013
For the Minister and by authorization
Ch. Tien-hsin
Secretary
Bureau of Consular Affairs
費用 fee: 400 NT\$400
G.35 978 5

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
Examination Lamp Mobile Stand Type EL-102, UL-102, UL-102S

Operating and Maintenance Instructions

 BenQ Medical Technology Corporation
7th Floor, No. 46, Zhou-Z Street, Nei-Hu,
Taipei 114, Taiwan
Customer Service e-mail :
Service@BenQMedicalTech.com
BenQMedicalTech.com
Customer Service Hotline :
Taiwan : 0800-015-533
China : +86 21-6327-7161~3 Ext. 812
International : +886 2-2601-7227 Ext. 5932

This instruction is intended to provide users with the proper information required for correct operation and maintenance of BenQ examination lamps. In order to avoid danger to patient or users and/or damage to the lamps, please review the contents of this instruction thoroughly and carefully before use.
This information contained in this instruction is subject to change without prior notice.

BenQ STANDARD LIMITED WARRANTY

BenQ warrants all new examination lamps sold by it directly or through a dealer or other authorized representatives, with expectation to replacement parts, spares, pads, and accessories, to be free from defects in material or workmanship, under normal use and service, for a period of one year. The warranty period shall begin with one year after shipment from BenQ. BenQ's responsibility and liability shall be limited to the replacement of any parts which we, BenQ, determine to be defective within the applicable warranty period, labor will not be included in the warranty coverage. No dealer, agent, or employee is authorized to extend or enlarge this warranty.

EC REP


ARTMED Ltd.

H-1138 Budapest, Viza u 7/c, B303, Hungary.

TABLE OF CONTENTS

1. Installation Instruction	Page 2
2. Operation Instruction	Page 3
3. Warning	Page 4
4. Cautions	Page 5
5. Maintenance Instruction	Page 6
6. Circuit Diagram	Page 7
7. Technical Data	Page 8

1. INSTALLATION INSTRUCTION

The BenQ Medical Technology Corp. model EL-102, UL-102 and UL-102S are facility examination lamp with mobile standing mount. There is no explosive gas mixture present. The unit is only for indoor use, and does not use it outdoors.

When you receive the well-packed operation lamp unit, please unpack the product immediately; in order to ensure the good quality of the BenQ Medical Technology Corp. product having been maintained during the transportation and the storage stages. After you are sure about the good quality of the product has been maintained well, you should check the packing list first. Check the content of the packing list be consistent with what you can find in the package. If there is any item missing or broken, please do not feel hesitate to contact with us. BenQ Medical Technology Corp., or the local distributor around your area for requiring or exchanging a new item.

Because all of the components of the operation lamp unit are already well assemble on the unit, just install as following figure.

2. OPERATION INSTRUCTION

First, the BenQ Medical Technology Corp. staff would like to thank you for choosing this examination lamp for your professional use. We hope the users can appreciate and get benefit from the high performance of the unit. Please note the following things:

- ◎ Read this operation instruction first before using the operation lamp, in order to understand the proper method of operation.
- ◎ Please save this operation instruction properly, and the user can refer to the content of the operation instruction when the user has doubt about products and does not know how to operate the unit.
- ◎ Whenever the users need technical assistance, please feel free to contact the local distributor of BenQ Medical Technology Corp., or directly contact the staff from the BenQ Medical Technology Corp.
- ◎ Emergency processing method: whenever the smoke or odor is found around the unit, please turn the lamp off immediately, and inform BenQ Medical Technology Corp. service personnel of dealing with this fault.
- ◎ This lamp unit is intended for use in health care facilities, such as hospital, clinics, and dental clinics, and laboratories, not for use in home-care facilities.

3. WARNING

The MUST about the safety operation of the examination lamp.

 Warning	
Use power with right rating.	Please make sure the tightness of the cover and connection is always in good condition.
Do not wear power cord. (any squeezing, twisting, wrapping, and pulling will induce breakage, and the hazard of fire and electric shock)	Do not use broken power cord, plug, or loosen socket. (may lead to hazard of electric shock, short-circuit, and hazard of fire)
Please not reform the unit.	Please not immerse the unit in water or rinse it with water. (may lead to hazard of electric shock)
Please clean operation lamp unit frequently.	Do not stare at the lamp when lamp is on.
Turn off / unplug and allow to cool before replacing lamp.	Lamp gets hot quickly, contact only switch / plug when turning on.
Do not touch hot lens 、 guard or enclosure.	Do not remain in light if skin feels warm.
Use only with a 35watt or smaller lamp.	Keep lamp away from materials that may burn.
Do not operate the lamp with a missing or damage guard 、 lamp containment barrier 、 or UV light filter.★	Do not touch the lamp at any time. Use a soft cloth. Oil from skin may damage lamp.

★ Remarks:

- Guard -- guard is a portion of portable lighting unit that prevents inadvertent contact with the lamp. May be integral with the UV filter or a portion of enclosure or shade.
- Lamp containment barrier -- lamp containment barrier is a part of a lighting unit that encloses the lamp.
- UV filter -- UV filter is a portion of the lighting unit that limits ultraviolet (UV) emission.

4. CAUTIONS

The NOTES when operating the examination lamp unit.

CAUTION	SOLUTION
Electric shock may occur while the lamp is being maintained or repaired.	Disconnect the power connection before further main tenance or repair.
Emission of high temperature from bulb may cause burnt hazard.	Do not change bulb when the cover of lamp-head is still warm. Wear the insulation glove to change the bulb until the bulb is cool.
Heat and high temperature emitted from lamp may cause fire. Do not place things easily on fire near the bulb to cause fire hazard.	Keep away from the lamp head as possible as you can.
Restriction pole keep body from lamp head is inclined to cause punch hurt.	Keep away from the lamp head.
Always hold the handle bar while the lamp head is being moved.	Always hold the handle bar when moving the lamp head.

5. MAINTENANCE INSTRUCTION

Trouble Shooting

SYMPTOMS	POSSIBLE CAUSES	RESOLUTIONS
A. The bulb go off.	(1) Lifespan of the bulb ends. (2) Broken wire. (3) Broken transformer. (4) Power cord is being cut. (5) Broken switch. (6) Broken fuse. (7) The fuse holder is loosen. (8) The current is wrong.	(1) Replace with spare bulbs.★ (2) Frequent occurrence can be found at connection. Check and clamp the loosen connection. (3) Inspect power supply and replace one. (4) Connect the broken wire. This may happen at arms during swing motion. (5) Replace new switch. (6) Check with short circuit if any and ensure the wiring is correct. Replace with new fuse. (7) Replace with new spring if the elasticity of old one is gone. (8) Examine power source, fuse, transformer and bulb.


★ Replacement of the bulb:

1. This light bulb socket shall be fitted with 12V 35W 12° .
2. It is recommended to replacement bulbs offered by manufacturers. The intensity, color temperature, beam angle ...etc of the bulb will affect the function.
3. Care should be taken while the bulb is being replaced. The bulb may be extremely hot.


6. CIRCUIT DIAGRAM

The following is enclosed copies of circuit diagram.

© DE-1102-EE01 (for UL-102)

© DE-1102-EE02 (for EL-102, UL-102S)

TECHNICAL DATA

Model	EL-102	UL-102	UL-102S
Item	Specification		
Number of Light Sources	1		
Size of Reflector (mm)	Ø 50		
Size of Lamphead (mm)	Ø 75		
Power (W)	35		
Intensity Control	NON-STEP		
Max. Intensity at 1m Distance (Lux)	13,000		
Color Temperature (°K)	3,500		
Control Method	MANUAL		
Minimum Ceiling Height	NA		
Weight (Kg)	5		
Voltage	240/230/220 /110/100V	120V	240/230/220 /110/100V

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