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Signed: [Signature]
Zaman Khan,
Associate Director Regulatory
Affairs EMEA and Pakistan,
Abbott GmbH & Co. KG.

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Toxo IgM Reagent Kit

en
Toxo IgM
07P47
G71285R06
B7P470

Read Highlighted Changes: Revised June 2019.

REF 07P4722

REF 07P4732

Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from these instructions.

NAME

Alinity i Toxo IgM Reagent Kit

INTENDED USE

The Alinity i Toxo IgM assay is a chemiluminescent microparticle immunoassay (CMIA) used for the qualitative detection of IgM antibodies to *Toxoplasma gondii* in human serum and plasma on the Alinity i analyzer.

SUMMARY AND EXPLANATION OF THE TEST

Toxoplasma gondii is an obligate intracellular protozoan parasite that infects most species of warm-blooded animals, including humans.¹ Toxoplasmosis is primarily acquired by ingestion of undercooked, infected meat; via oocysts from fecally contaminated hands, food and water; and maternally through transplacental transmission.² In addition, transmission associated with organ transplantation and during blood transfusion has been reported, although the risk of transmission through blood transfusion is extremely low.³

Acquired infection with *Toxoplasma gondii* in healthy individuals is commonly asymptomatic, however 10-20% of patients with acute infection may develop lymphadenopathy.⁴

Severe infections can occur in AIDS patients and adults immunocompromised by cancer chemotherapy or transplant recipients receiving immunosuppressive treatment. These infections can be fatal. Toxoplasmic encephalitis is the most common presentation and is the most frequent cause of focal central nervous system lesions in AIDS patients.⁵

Primary infection during pregnancy can result in transplacental transmission of the parasite resulting in congenital infection. The risk of congenital infection is lowest (10-25%) if acute maternal infection occurs during the first trimester and highest (60-90%) if it occurs during the third trimester.²

Severity of congenital infection is greatest when maternal infection is acquired during early pregnancy. Common outcomes of congenital toxoplasmosis include chorioretinitis, intracranial calcifications, and hydrocephalus. The majority of infants infected later in pregnancy are asymptomatic at birth, with sequelae occurring later in life.

Early treatment after prenatal diagnosis of *Toxoplasma gondii* infection has been shown to reduce the frequency and severity of congenital toxoplasmosis.⁶ Serological tests can be used to identify seronegative women, who then should be monitored during pregnancy.

The presence of IgG antibodies to *Toxoplasma gondii* indicates that infection has occurred but does not distinguish between recent and past infection. IgM antibodies are detected in individuals with a recently acquired infection, but antibodies may persist for up to 18 months post-infection.²

To differentiate between a recently acquired and past infection, IgM and IgG positive specimens should be tested for IgG avidity. A high avidity index for IgG antibodies is a strong indication that an infection took place more than 4 months ago. Low avidity results cannot be used to diagnose an acute toxoplasmosis.

Toxo IgG	Toxo IgM	Toxo IgG Avidity	May indicate... / Testing recommendation
nonreactive	nonreactive	N/A	no infection
nonreactive	reactive	N/A	obtain new sample 2-3 weeks after initial sample and test for Toxo IgG and Toxo IgM
reactive	nonreactive	high avidity	past infection. Strong indication that an infection took place more than 4 months ago
reactive	reactive	low avidity	obtain new sample 3 weeks after initial sample and test for Toxo IgG and Toxo IgM
reactive	reactive	high avidity	past infection. Strong indication that an infection took place more than 4 months ago

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

This assay is a two-step immunoassay for the qualitative detection of IgM antibodies to *Toxoplasma gondii* in human serum and plasma using chemiluminescent microparticle immunoassay (CMIA) technology.

Pre-diluted sample and anti-human IgM mouse monoclonal antibody coated paramagnetic microparticles are combined and incubated. Together with IgM antibodies of other specificities, anti-Toxo specific IgM present in the sample binds to the anti-human IgM mouse monoclonal antibody coated microparticles, forming an antibody-antibody complex. The mixture is washed. A conjugate complex consisting of an acridinium-labeled anti-Toxo p30 antigen mouse monoclonal F(ab')₂ fragment and native *Toxoplasma gondii* lysate, containing the p30 antigen, is added to create a reaction mixture and incubated. This conjugate complex is bound by anti-Toxo specific IgM that has been captured by the anti-human IgM mouse monoclonal antibody coated microparticles, forming an antibody-antibody-conjugate complex. Following a wash cycle, Pre-Trigger and Trigger Solutions are added.

The resulting chemiluminescent reaction is measured as relative light units (RLUs). There is a direct relationship between the amount of anti-Toxo IgM in the sample and the RLUs detected by the system optics.

The presence or absence of anti-Toxo IgM in the sample is determined by comparing the chemiluminescent RLU in the reaction to the cutoff RLU determined from an active calibration.

For additional information on system and assay technology, refer to the Alinity ci-series Operations Manual, Section 3.

REAGENTS

Kit Contents

Alinity i Toxo IgM Reagent Kit 07P47

NOTE: Some kit sizes are not available in all countries. Please contact your local distributor.

Volumes (mL) listed in the table below indicate the volume per cartridge.

REF	07P4722	07P4732
Tests per cartridge	100	500
Number of cartridges	2	2
Tests per kit	200	1000
MICROPARTICLES	6.6 mL	27.0 mL
CONJUGATE	6.1 mL	26.5 mL

MICROPARTICLES Anti-human IgM (mouse, monoclonal) antibody coated microparticles in TRIS buffer with protein stabilizers (bovine and goat) and detergent. Minimum concentration: 0.08% solids. Preservatives: antimicrobial agents.

CONJUGATE Conjugate complex consisting of acridinium-labeled anti-Toxoplasma p30 antigen (mouse, monoclonal) antibody and native *Toxoplasma gondii* lysate in phosphate buffer with protein stabilizers (bovine) and detergent. Minimum concentration: 25 µg/mL. Preservative: sodium azide.

Warnings and Precautions

- **IVD**
- For *In Vitro* Diagnostic Use

Safety Precautions

CAUTION: This product contains human-sourced and/or potentially infectious components. Refer to the **REAGENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.⁷⁻¹⁰

The following warnings and precautions apply to: CONJUGATE	
Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

Reagent Handling

- Upon receipt, gently invert the unopened reagent kit by rotating it over and back for a full 180 degrees, 5 times with green label stripe facing up and then 5 times with green label stripe facing down. This ensures that liquid covers all sides of the bottles within the cartridges. During reagent shipment, microparticles can settle on the reagent septum.
 - Place a check in the square on the reagent kit to indicate to others that the inversions have been completed.
- After mixing, place reagent cartridges in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- If a reagent cartridge is dropped, place in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- Reagents are susceptible to the formation of foam and bubbles. Bubbles may interfere with the detection of the reagent level in the cartridge and cause insufficient reagent aspiration that may adversely affect results.

For a detailed discussion of reagent handling precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 7.

Reagent Storage

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	Store in upright position. If cartridge does not remain upright, gently invert the cartridge 10 times and place in an upright position for 1 hour before use.
Onboard	System Temperature	30 days	
Opened	2 to 8°C	Until expiration date	Store in upright position. If cartridge does not remain upright during storage, discard the cartridge. Do not reuse original reagent caps or replacement caps due to the risk of contamination and potential to compromise reagent performance.

Reagents may be stored on or off the system. If removed from the system, store reagents with new replacement caps in an upright position at 2 to 8°C. For reagents stored off the system, it is recommended that they be stored in their original trays or boxes to ensure they remain upright.

For information on unloading reagents, refer to the Alinity ci-series Operations Manual, Section 5.

Indications of Reagent Deterioration

Deterioration of the reagents may be indicated when a calibration error occurs or a control value is out of the specified range. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary.

For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.

INSTRUMENT PROCEDURE

The Alinity i Toxo IgM assay file must be installed on the Alinity i analyzer prior to performing the assay.

The Alinity ci-series software version 2.0 or higher must be installed prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the Alinity ci-series Operations Manual, Section 2.

For information on printing assay parameters, refer to the Alinity ci-series Operations Manual, Section 5.

For a detailed description of system procedures, refer to the Alinity ci-series Operations Manual.

Alternate Result Units

Edit assay parameter "Result Units" to select an alternate unit.

Conversion formula:

(Default result unit) / (Conversion divisor) = (Alternate result unit)

Default Result Unit	Conversion Divisor	Alternate Result Unit
Index	0.60	S/CO

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Types

The specimen types listed below were verified for use with this assay on the ARCHITECT i System.

Other specimen types and collection tube types have not been verified with this assay.

Specimen Types	Collection Tubes
Serum	Serum
	Serum separator
Plasma	Plasma separator (lithium heparin)
	Potassium EDTA
	Sodium heparin
	Lithium heparin

- Performance has not been established for the use of cadaveric specimens or the use of bodily fluids other than human serum and plasma.
- Liquid anticoagulants may have a dilution effect resulting in lower concentration values for individual specimens.
- Sodium citrate tubes cannot be used with the Alinity i Toxo IgM assay.
- The instrument does not provide the capability to verify specimen types. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

Specimen Conditions

- Do not use:
 - heat-inactivated specimens
 - pooled specimens
 - grossly hemolyzed specimens
 - specimens with obvious microbial contamination
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

Preparation for Analysis

- Follow the tube manufacturer's processing instructions for collection tubes. Gravity separation is not sufficient for specimen preparation.
- Specimens should be free of bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

To ensure consistency in results, recentrifuge specimens prior to testing if

- they contain fibrin, red blood cells, or other particulate matter
- they require repeat testing.

NOTE: If fibrin, red blood cells, or other particulate matter are observed, mix by low speed vortex or by inverting 10 times prior to recentrifugation.

Prepare frozen specimens as follows:

- Frozen specimens must be completely thawed before mixing.
- Mix thawed specimens thoroughly by low speed vortex or by inverting 10 times.
- Visually inspect the specimens. If layering or stratification is observed, mix until specimens are visibly homogeneous.
- If specimens are not mixed thoroughly, inconsistent results may be obtained.
- Recentrifuge specimens.

Recentrifugation of Specimens

- Transfer specimens to a centrifuge tube and centrifuge at a minimum of 100 000 g-minutes.
 - Examples of acceptable time and force ranges that meet this criterion are listed in the table below.
- Centrifugation time using alternate RCF values can be calculated using the following formula:

$$\text{Minimum Centrifugation time (minutes)} = \frac{100\,000 \text{ g-minutes}}{\text{RCF}}$$

Recentrifugation Time (Minutes)	RCF (x g)	g-Minutes
10	10 000	100 000
20	5000	100 000
40	2500	100 000

$$\text{RCF} = 1.12 \times r_{\text{max}} (\text{rpm}/1000)^2$$

RCF - The relative centrifugal force generated during centrifugation.

rpm - The revolutions per minute of the rotor on which the specimens are being spun (usually the digital readout on the centrifuge will indicate the rpm).

Centrifugation Time - The time should be measured from the time the rotor reaches the required RCF or rpm to the time it begins decelerating.

r_{max} - Radius of the rotor in millimeters. NOTE: If custom tube adapters (i.e., adapters not defined by the centrifuge manufacturer) are used, then the radius (r_{max}) should be manually measured in millimeters and the RCF calculated.

g-minutes - The unit of measure for the product of RCF (x g) and centrifugation time (minutes).

- Transfer clarified specimen to a sample cup or secondary tube for testing. For centrifuged specimens with a lipid layer, transfer only the clarified specimen and not the lipemic material.

Specimen Storage

Specimen storage conditions were verified on the ARCHITECT i System.

Specimen Type	Temperature	Maximum Storage Time	Special Instructions
Serum/ Plasma	15 to 30°C	3 days	Specimens may be stored on or off the clot, red blood cells, or separator gel.
	2 to 8°C	14 days	Specimens may be stored on or off the clot, red blood cells, or separator gel.

Remove serum or plasma from the clot, red blood cells, or separator gel if stored longer than the maximum 15-30°C or 2-8°C storage time and store frozen at -20°C or colder.

Avoid more than 6 freeze/thaw cycles.

Specimen Shipping

Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.

PROCEDURE

Materials Provided

07P47 Alinity i Two Igm Reagent Kit

Materials Required but not Provided

- Alinity i Two Igm assay file
- 07P4701 Alinity i Two Igm Calibrator
- 07P4710 Alinity i Two Igm Controls or other control material
- Alinity Trigger Solution
- Alinity Pre-Trigger Solution
- Alinity i-series Concentrated Wash Buffer

For information on materials required for operation of the instrument, refer to the Alinity i-series Operations Manual, Section 1.

For information on materials required for maintenance procedures, refer to the Alinity i-series Operations Manual, Section 9.

Assay Procedure

For a detailed description of how to run an assay, refer to the Alinity i-series Operations Manual, Section 5.

- If using primary or aliquot tubes, refer to the Alinity i-series Operations Manual, Section 4 to ensure sufficient specimen is present.
- To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
- Maximum number of replicates sampled from the same sample cup: 10
 - Priority:
 - Sample volume for first test: 70 µL
 - Sample volume for each additional test from same sample cup: 20 µL
 - ≤ 3 hours on the reagent and sample manager:
 - Sample volume for first test: 150 µL
 - Sample volume for each additional test from same sample cup: 20 µL
 - > 3 hours on the reagent and sample manager:
 - Replace with a fresh aliquot of sample.
- Refer to the Alinity i Two Igm calibrator package insert and/ or Alinity i Two Igm control package insert for preparation and usage.
- For general operating procedures, refer to the Alinity i-series Operations Manual, Section 5.
- For optimal performance, it is important to perform routine maintenance as described in the Alinity i-series Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Sample Dilution Procedures

Samples cannot be diluted for the Alinity i Two Igm assay.

Calibration

For instructions on performing a calibration, refer to the Alinity i-series Operations Manual, Section 5.

Each assay control must be tested to evaluate the assay calibration. Once a calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:

- A reagent kit with a new lot number is used.
- Daily quality control results are outside of statistically-based quality control limits used to monitor and control system performance, as described in the Quality Control Procedures section of this package insert.
 - If statistically-based quality control limits are not available, then the calibration should not exceed a 30-day limit for recalibration frequency.

This assay may require recalibration after maintenance to critical parts or subsystems or after service procedures have been performed.

Quality Control Procedures

The recommended control requirement for the Alinity i Two Igm assay is that a single sample of each control level be tested once every 24 hours each day of use.

Additional controls may be tested in accordance with local, state, and/or federal regulations or accreditation requirements and your laboratory's quality control policy.

To establish statistically-based control limits, each laboratory should establish its own concentration target and ranges for new control lots at each clinically relevant control level. This can be accomplished by assaying a minimum of 20 replicates over several (3-5) days and using the reported results to establish the expected average (target) and variability about this average (range) for the laboratory. Sources of variation that should be included in this study in order to be representative of future system performance include:

- Multiple stored calibrations
- Multiple reagent lots
- Multiple calibrator lots
- Multiple processing modules (if applicable)
- Data points collected at different times of the day

Refer to published guidelines for information or general control recommendation, for example Clinical and Laboratory Standards Institute (CLSI) Document C24-A5 or other published guidelines, for general quality control recommendations.¹¹

- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary. For troubleshooting information, refer to the Alinity i-series Operations Manual, Section 10.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

Quality Control Guidance

Refer to "Basic QC Practices" by James O Westgard, Ph.D. for guidance on laboratory quality control practices.¹²

Verification of Assay Claims

For protocols to verify package insert claims, refer to Verification of Assay Claims in the Alinity i-series Operations Manual.

RESULTS

Calculation

The Alinity i analyzer calculates results for the Alinity i Two Igm assay using the ratio of the sample RLU to the Calibrator 1 Mean RLUI (Index) for each specimen and control.

The Calibrator 1 Mean RLUI is stored for each reagent lot calibration.
Index = Sample RLUI/Calibrator 1 Mean RLUI

Interpretation of Results

The cutoff is 0.60 Index (1.60 S/CO).

Initial Results

Specimens with results	Instrument Interpretation	Repeat Procedure
< 0.60 Index ($< 0.60 \text{ S/CO}$)	Nonreactive for IgM antibodies to Hepatitis virus	No repeat required.
$0.60 \leq X < 1.60$ Index ($0.60 \leq X < 1.60 \text{ S/CO}$)	Grayzone	It is recommended to take a second sample within a reasonable period of time (e.g. two weeks) and repeat Alinity i Two Igm testing.
≥ 0.60 Index ($\geq 1.60 \text{ S/CO}$)	Reactive for IgM antibodies to Hepatitis virus	No repeat required.

Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the Alinity ci-series Operations Manual, Section 5.

LIMITATIONS OF THE PROCEDURE

- If the Alinity i Toxo IgM results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.
- For diagnostic purposes, results should be used in conjunction with other data; e.g., results of other tests (Toxo IgG, Toxo IgG Avidity), clinical impressions, etc.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays. Patients routinely exposed to animals or to animal serum products can be prone to this interference, and anomalous values may be observed. Additional information may be required for diagnosis.¹³
- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits such as Alinity i Toxo IgM that employ mouse monoclonal antibodies. Additional information may be required for diagnosis.^{14, 15}
- Specimens from patients with high levels of IgM, e.g. from patients with multiple myeloma, may show depressed values when tested with μ -capture format assays.
- Human plasma collected in sodium citrate tubes cannot be used with this assay, since the Alinity i Toxo IgM results may be affected by this tube type.

SPECIFIC PERFORMANCE CHARACTERISTICS

Representative performance data are provided in this section. Results obtained in individual laboratories may vary. The Alinity i analyzer and the ARCHITECT i System utilize the same reagents and sample/reagent ratios. Unless otherwise specified, all studies were performed on the Alinity i analyzer.

Precision

Within-Laboratory Precision

A study was performed based on guidance from CLSI EP05-A2.¹⁶ Testing was conducted using 3 lots of the Alinity i Toxo IgM Reagent Kit, 3 lots of the Alinity i Toxo IgM Calibrator, and 3 lots of the Alinity i Toxo IgM Controls and 1 instrument. Two controls and 2 human plasma panels were assayed in a minimum of 2 replicates at 2 separate times per day on 20 different days.

Sample	n	Within-Run (Repeatability)		Within-Laboratory (Total) ^a	
		Mean S/CO	SD	SD (Range ^b)	%CV (Range ^b)
Negative Control	357	0.10	0.006	0.011 (0.009 - 0.014)	10.4 (9.6 - 11.0)
Positive Control	359	1.48	0.039	0.048 (0.044 - 0.052)	3.2 (3.0 - 3.5)
Panel 1	358	0.38	0.012	0.014 (0.013 - 0.016)	3.7 (3.5 - 3.9)
Panel 2	360	0.79	0.025	0.028 (0.025 - 0.030)	3.5 (3.3 - 3.7)

^a Includes within-run, between-run, and between-day variability.
^b Maximum and minimum SD or %CV for each reagent lot and instrument combination.

Sample	n	Within-Run (Repeatability)		Within-Laboratory (Total) ^a	
		Mean S/CO	SD	SD (Range ^b)	%CV (Range ^b)
Negative Control	357	0.17	0.010	0.019 (0.016 - 0.023)	11.0 (10.6 - 11.1)
Positive Control	359	2.46	0.065	0.079 (0.073 - 0.087)	3.2 (3.0 - 3.5)
Panel 1	358	0.64	0.020	0.024 (0.022 - 0.027)	3.8 (3.6 - 3.9)
Panel 2	360	1.31	0.041	0.046 (0.042 - 0.050)	3.5 (3.3 - 3.7)

^a Includes within-run, between-run, and between-day variability.
^b Maximum and minimum SD or %CV for each reagent lot and instrument combination.

Resolved Relative Specificity

This study was performed on the ARCHITECT i System. The ARCHITECT Toxo IgM assay is designed to have a resolved relative specificity comparable to a comparator Toxo IgM assay. A study was performed at one internal and one external evaluation site. From the 2772^{***} specimens evaluated to assess resolved relative specificity 36 specimens were concordant reactive and additional three were confirmed positive after discordant resolution and therefore excluded from the specificity calculation.
^{***} Note: Specimens that could not be resolved or showed grayzone result interpretation on any assay being compared or used for discordant resolution were not included in the evaluation of resolved relative specificity.

Data from this study are summarized in the following table.

Resolved Relative Specificity

Sample Type	ARCHITECT Toxo IgM		Comparator Toxo IgM Assay	
	Observed	Lower 95% Confidence Limit	Observed	Lower 95% Confidence Limit
Pregnant Women	99.95% (1987/1988)	99.72%	99.95% (1987/1988)	99.72%
Diagnostic / Hospital Patients	100% (451/451)	99.19%	100% (451/451)	99.19%
Blood Donors (Serum)	100% (154/154)	97.63%	100% (154/154)	97.63%
Blood Donors (Plasma)	98.57% (138/140)	94.93%	100% (140/140)	97.40%
Total	99.89% (2730/2733)	99.68%	99.96% (2732/2733)	99.80%

Negative Percent Agreement

A study was performed to determine the negative percent agreement between the Alinity i Toxo IgM assay and a commercially-available Toxo IgM assay by testing 1028 specimens on a population consisting of pregnant women, hospitalized individuals, and blood donors. The specimens were tested in singlicate using 3 Alinity i Toxo IgM Reagent Kit lots, 2 lots of Alinity i Toxo IgM Calibrator, and 1 lot of Alinity i Toxo IgM Controls. The specimens were also tested in singlicate using a commercially-available Toxo IgM assay. Specimens that showed a grayzone result interpretation on either the Alinity i Toxo IgM assay or the commercially-available Toxo IgM assay were not included in the evaluation. The negative percent agreement was 100.00% (i.e., 1028/1028, 1027/1027, and 1026/1026) on all 3 reagent lots tested. The lower 95% confidence interval (of the assumed comparator assay specificity of 100%) was 99.64%.

Seroconversion Sensitivity

This study was performed on the ARCHITECT i System.

The ARCHITECT T_{oxo} IgM assay is designed to show seroconversion sensitivity comparable to a comparator T_{oxo} IgM assay. A total of 122 bleeds from 59 different individuals seroconverting during acute toxoplasmosis infection were tested and demonstrated acceptable seroconversion detection. Data from seven representative seroconverting individuals are shown in the following table.

Sample ID	Months after last negative bleed	ARCHITECT T _{oxo} IgM [Index]	Comparator T _{oxo} IgM [Index]	WADA (T _{oxo} -M) [Index]	Subtle-Pathman IgM Test [S/Co] (S/Co)	MS Application Test (S/Co)
Group 1	0.00 - 0.50	0.000 - 0.500	0 - 5	ANA	1	
Positive cutoff	0.00	0.000	0	2	2	
30044001	0.0	0.06	0.061	0	< 2	< 1
30044002	1.2	0.66	0.361	4	5	7
30044003	2.1	0.80	0.600	10	800	94
30044004	2.2	0.77	0.542	10	800	94
30044005	4.3	0.30	0.193	1	1000	200
30044016	6.8	0.84	0.800	0	< 2	< 1
30044017	1.4	1.10	1.000	11	30	2
30044018	1.8	1.84	0.907	10	30	4
30044019	4.1	0.62	0.300	0	10	5
30044023	0.0	0.04	0.210	0	< 2	< 1
30044034	2.6	1.90	2.000	12	400	64
30044035	7.5	0.05	0.067	0	800	100
30044073	0.0	0.04	0.076	0	< 2	< 1
30044074	0.9	1.32	1.125	12	5	7
30044075	1.4	2.26	1.730	12	200	10
30044076	3.0	1.12	0.877	12	100	5
30044086	0.0	0.20	0.407	7	2	< 1
30044087	0.5	10.30	7.074	12	40	5
30044088	1.3	0.25	0.404	12	400	120
30044089	2.3	0.95	0.300	12	400	90
30044090	0.0	0.05	0.061	0	< 2	< 1
30044091	1.2	0.95	0.105	12	30	2
30044092	1.5	0.72	0.679	12	200	10
30044093	4.7	2.00	1.700	12	400	90
30044119	6.0	0.65	0.115	0	< 2	< 1
30044119	1.0	0.76	0.704	12	2	7
30044120	1.8	0.80	0.536	12	200	64
30044121	2.5	0.80	1.010	12	1000	100

This study was performed on the Alinity i analyzer.

To determine the seroconversion sensitivity, 1 seroconversion panel obtained from a commercial vendor was tested on the Alinity i analyzer using the Alinity i T_{oxo} IgM assay. The panel results were evaluated against a commercially-available T_{oxo} IgM assay and data are summarized in the following table.

Panel ID	Days Since Test Bleed	Alinity i T _{oxo} IgM (Index)			Commercially-Available T _{oxo} IgM Assay (Index)
		Reagent Lot 1	Reagent Lot 2	Reagent Lot 3	
451657152	0	0.20	0.20	0.20	0.22
	7	1.04	1.20	1.11	1.00
	9	2.22	2.00	2.00	2.00
	14	0.82	7.01	5.80	0.76
	16	0.70	7.50	5.70	0.60
	23	0.27	0.80	0.84	0.42
	27	4.00	0.21	0.15	4.00
	30	4.40	0.70	4.74	4.00
	40	3.40	4.40	3.34	3.00
	63	3.11	4.00	3.27	2.00
	66	2.40	3.20	2.80	2.00
	86	2.10	2.74	2.40	2.27
	90	1.70	2.40	2.10	1.00
	95	1.71	2.01	1.80	1.00
	100	1.80	2.47	1.80	1.70
	104	1.80	2.21	1.80	1.77
	107	1.84	2.12	1.80	1.80
	112	1.70	1.90	1.70	1.82
	114	1.40	1.90	1.82	1.44
	118	1.44	2.01	1.82	1.43

Other Specimen Conditions or Disease States

This study was performed on the ARCHITECT i System.

Studies were performed to evaluate other potential interfering disease states on the ARCHITECT T_{oxo} IgM assay. Sample categories were tested both unspiked and spiked with anti-T_{oxo} IgM positive plasma.

A total of 167 unspiked specimens and a total of 165 spiked anti-T_{oxo} IgM reactive specimens were tested from the following categories: Anti-nuclear antibodies (ANA), anti-dsDNA antibodies, Rheumatoid Factor, Herpes Simplex Virus 1 (anti-HSV-1 positive), Herpes Simplex Virus 2 (anti-HSV-2 positive), Epstein-Barr Virus (anti-EBV positive), Measles, Parvovirus B19 (anti-B19 virus positive), Varicella Zoster Virus (anti-VZV positive), Rubella Virus (anti-Rubella positive), Cytomegalovirus (anti-CMV positive), Hyperpolyclonal IgG, Hyperpolyclonal IgM, Monoclonal IgG, Monoclonal IgM, Human anti-mouse antibodies (HAMA), Influenza vaccine recipients, and Syphilis. ARCHITECT T_{oxo} IgM showed expected qualitative results in all categories with the exception of the categories summarized in the following table.

Category	N tested	ARCHITECT T _{oxo} IgM		
		Nonreactive	Reactive	Group 1
Anti-nuclear Antibodies (ANA) Unspiked	10	8	1	1
Hyperpolyclonal IgM Spiked Reactive	10	1	9	1
Monoclonal IgM Spiked Reactive	5	4	0	1

After discordant resolution the unspiked ANA specimen tested reactive on ARCHITECT T_{oxo} IgM could not be resolved. Specimens from patients with high levels of IgM, e.g. from patients with multiple myeloma, may show depressed values when tested with μ -capture format assays. Refer to the LIMITATIONS OF THE PROCEDURE section of this package insert.

Interference

This study was performed on the ARCHITECT i System.

Potentially Interfering Endogenous Substances

No interference was observed between experimental controls and nonreactive or reactive specimens tested with elevated levels of compounds listed in the table below.

Potentially Interfering Substance	Interferent Level
Bilirubin	≤ 20 mg/dL
Triglycerides	≤ 3000 mg/dL
Protein	≤ 12 g/dL
Hemoglobin	≤ 500 mg/dL

BIBLIOGRAPHY

- Willis MS, Southern P, Latimer MJ. Toxoplasma infection. Making the best use of laboratory tests. *Infect Med* 2002; 19: 522-532.
- Jones JL et al. Congenital Toxoplasmosis: A Review. *CME Review Article* Vol. 56, Number 5 2001; 296-305.
- Shulman IA. Parasitic infections and their impact on blood donor selection and testing. *Arch Pathol Lab Med*. 1994; 118: 366-370.
- Montoya J.G. Laboratory Diagnosis of Toxoplasma gondii Infection and Toxoplasmosis. *Journal of Infectious Diseases* 2002; 185 (Suppl 1) 73-82.
- Israelski DM, Remington JS. Toxoplasmosis in the non-AIDS immunocompromised host. *Curr Clin Top Infect Dis* 1993; 13: 322-356.
- Wong SY and Remington. Toxoplasmosis in Pregnancy. *Clinical Infectious Diseases* 1994; 18:853-862.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Clinical and Laboratory Standards Institute (CLSI). *Statistical Quality Control for Quantitative Measurement Procedures: Principles and Definitions; Approved Guideline—Third Edition*. CLSI Document C24-A3. Wayne, PA: CLSI; 2006.
- Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
- Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261-264.
- Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879-885.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
	Conjugate
	Contains Sodium Azide. Contact with acids liberates very toxic gas.
	Inversions Performed
	In Vitro Diagnostic Medical Device
	Batch code/Lot number
	Microparticles
	Product of Germany
	List Number
	Serial number

Alinity and ARCHITECT are trademarks of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

0123

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Revised June 2019.
©2016, 2019 Abbott Laboratories

ADDENDUM

to the instruction for Use Alinity I Toxo IgM Reagent Kit for Russia

SHIPPING AND STORAGE CONDITIONS

Recommended shipping conditions for devices are 2-8°C. Medical devices must be stored in refrigerators/refrigerator chambers to ensure regulated temperature conditions with a daily temperature record.

Medical devices must be carried out in single-use thermal containers containing coolant elements, or in reusable thermo containers with automatically maintained temperature or in refrigerated trucks using thermal indicators.

Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an any kind of enclosed carrier vehicle in accordance with shipping rules, applicable to the respective type of transport.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as intended. Manufacturer set the maximum shelf life of

Alinity I Toxo IgM Reagent Kit – 12 months from date of manufacturing.

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste". Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

98/78/EC, 1272/2008/EC (CLP), 1907/2006/EC (REACH), 1069/2009, Regulation (EC) 142/2011, CLASS 2013, GHS, WHSR Part 7.1, EN ISO 9001, ISO 13485, EN ISO 14871, EN ISO 13641, EN ISO 13612, EN ISO 15223-1, EN 980 EN ISO 18113-1, EN ISO 18113-2, EN ISO 18113-3, EN ISO 23840, EN ISO 20776-1, EN 13675, EN ISO 17511

POTENTIAL USER

The Alinity I Toxo IgM Reagent Kit is intended to be used in a variety of clinical and laboratory conditions, including, without limitation, hospitals, diagnostic centers, reference laboratories and blood transfusion facilities around the world, whose personnel have been trained in the use of the Alinity I analyzer.

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device is intended for one-off usage

MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

INDICATIONS AND CONTRAINDICATIONS FOR USE

Not applicable.

TECHNICAL CHARACTERISTICS

- Microparticles: pH-value at 20 °C: 7,4
Density at 20 °C: 1.007 g/cm³
- Conjugate: pH-value at 20 °C: 6,3
Density at 20 °C: 1.014 g/cm³

COMPATIBILITY RESTRICTIONS OF MEDICAL DEVICE

No available data.

COLOUR AND ODOUR

Microparticles: Colourless and odourless
Conjugate: Colourless and odourless

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- Willis MS, Southern P, Latimer MJ. Toxoplasma infection. Making the best use of laboratory tests. *Infect Med* 2002; 19: 522-532.
- Jones JL et al. Congenital Toxoplasmosis: A Review. *CME Review Article* Vol. 56, Number 5 2001; 296-305.
- Shulman IA. Parasitic infections and their impact on blood donor selection and testing. *Arch Pathol Lab Med.* 1994; 118: 366-370.
- Montoya J.G. Laboratory Diagnosis of Toxoplasma gondii Infection and Toxoplasmosis. *Journal of Infectious Diseases* 2002; 185 (Suppl 1) 73-82.
- Israelski DM, Remington JS. Toxoplasmosis in the non-AIDS immunocompromised host. *Curr Clin Top Infect Dis* 1993; 13: 322-356.
- Wong SY and Remington. Toxoplasmosis in Pregnancy. *Clinical Infectious Diseases* 1994; 18:853-862.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Clinical and Laboratory Standards Institute (CLSI). *Statistical Quality Control for Quantitative Measurement Procedures: Principles and Definitions; Approved Guideline—Third Edition*. CLSI Document C24-A3. Wayne, PA: CLSI; 2006.
- Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
- Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261-264.
- Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879-885.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.

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