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Thank you for placing your trust in the name KARL STORZ. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high quality piece of equipment from KARL STORZ.

This instruction manual, in conjunction with the instruction manual for the camera control unit used, is intended to serve as an aid in the proper connection and operation of the camera head. All required details and all actions on your part are clearly presented and explained. Please read this manual carefully before proceeding with the equipment. Keep this instruction manual available for ready reference in a convenient location near the equipment.

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It is recommended to check the suitability of the product for the intended procedure prior to use.

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1 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

1.1 Explanation of warnings and cautions

Please read this manual and follow the instructions carefully. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system. To make these words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician is at risk.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

NOTE: A Note indicates special information about operating the device or clarifies important information.

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1.2 General safety information

-  **WARNING:** Read this instruction manual thoroughly before using the equipment.
-  **WARNING:** Check that the equipment is functioning properly prior to each surgical intervention. Ensure that the correct video image appears on all video monitors before beginning each intervention. It is recommended that a replacement system is kept available.
-  **WARNING:** Perceived electric shocks may be experienced when using the camera system in conjunction with electrosurgical units.
-  **WARNING:** When using high frequency surgical accessories, care must be taken to avert contact between the HF surgical device and the endoscope or endoscopic camera. Inadvertent contact may represent a thermal hazard to the surgeon.
-  **WARNING:** When energized endoscopes are used with energized endotherapy devices, patient leakage currents may be additive. It should be noted that this is particularly important if a Type CF applied part endoscope is used, in which case a Type CF applied part energized endotherapy device should be used in order to minimize total patient leakage current.
-  **WARNING:** Always ensure the connector is completely clean and free of debris. If moisture is present, dry thoroughly with sterile towel prior to insertion.

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WARNING: HF-generating equipment is dangerous and cannot be allowed to make contact with the camera head or system via the accessories or endoscopes connected to the camera head or system.

WARNING: Risk of blindness. Never look into the open end of a light cable.

WARNING: Risk of burns. Contact with the open end of a light cable connected to a light source may cause burns.

WARNING: Test all equipment prior to each surgical use. In the event that the image becomes unusable during surgery, the camera may be disengaged from the endoscope and the procedure continued optically. If this is not possible, it is up to the discretion of the surgeon how best to proceed. Availability of a spare system is recommended.

WARNING: High energy radiated light through endoscopes may give rise to high temperatures in front of the light outlet and to the tip of the endoscope. To minimize the risk of burns, always adjust the light source to the minimum illumination intensity necessary to achieve optimum illumination of the endoscopic scene when coupled to the video camera.

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2 Intended Use Camera Heads

KARL STORZ camera heads, in combination with the corresponding camera control unit, light source, monitor and telescopes or microscopes, are used for real-time display of the endoscopic or open operating field in diagnostic or surgical procedures.

KARL STORZ camera heads have no direct contact with the patient.

2.1 Contraindications

No contraindications relating directly to the medical device are currently known. The responsible physician must decide whether the foreseen application is admissible based on the general condition of the patient.

2.2 User qualification

The H1 camera head may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in the use of the H1 camera head.

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



- ① Instrument coupler
- ② Focusing ring/knob
- ③ User control button: for scrolling (downwards) through the menu or activating a camera function
- ④ User control button: for scrolling (upwards) through the menu or activating a camera function
- ⑤ User control button: for accessing the menu/ making a selection

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4 Symbols employed

	Fragile, handle with care
	Keep dry
	Temperature limit
	Serial number
	Number of products in the product packaging
	CE marking
	top
	Consult instructions for use
	Manufacturer
	Federal (USA) law restricts this device to sale by or on the order of a physician.

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5 Operating instructions for the H1 camera head

NOTE: While operating with surgical lasers intensive laser light may occur and affect image visualization.

5.1 Frequency of the H1 camera head (25/30 Hz)

The H1 camera head is suitable for operation at either 25 Hz or 30 Hz. The camera head must be selected to suit the local conditions.

NOTE: The switchover is described in the instruction manual of the respective camera control unit used.

5.2 Operating the H1 camera head

5.2.1 Connecting an endoscope and light cable

The H1 camera head has an integrated instrument coupler and integrated telescope. To connect an endoscope to the instrument coupler, turn the outer ring of the coupler clockwise and insert the endoscope eyepiece. To prevent rotation of the endoscope, the outer ring of the instrument coupler must be turned counterclockwise until tight. Connect the light cable to the endoscope eyepiece by tightening the knurled screw on the light cable through a quarter turn.

CAUTION: Heat from the light source can cause disinfectant residue on the light cable or on the light inlet of the endoscope to burn. This can have a detrimental effect on the light conduction and thus the image quality.

NOTE: For further information regarding light cables, please refer to the manual „Fiberoptic light cables“ (97000213), which can be downloaded or requested at www.karlstorz.com.

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5.2.2 Adjusting the focus

The focus on the H1 camera head can be adjusted by turning the focusing ring/knob ② on the camera lens.

5.2.3 Operation of camera head buttons

• In combination with TELE PACK +

Functions can be freely assigned to the control buttons (see operating instructions TELE PACK + 96286007DER).

• In combination with C-HUB® II or C-MAC® monitors

The following functions are permanently assigned to the camera head buttons and cannot be freely assigned.

Left control button: Image capture

Right control button: Start/stop video recording

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7 Reprocessing (not validated for US customers)

In the following chapter, you will find information on care and cleaning of the H1 camera head.

7.1 Care and handling

To maximize the life of your camera, please follow the guidelines given below for the correct handling and care of the camera head:

DO NOT PLUG THE CAMERA HEAD CABLE INTO THE PROCESSOR RECEPTACLE IF IT IS MOIST OR NOT CLEAN. Always ensure that the connector is completely clean. Remove any moisture present with a sterile cloth before connecting.

DO NOT PICK UP THE CAMERA HEAD BY ITS CABLE. Always pick up the camera head first, and then pick up and support the cable and connector. The cable contains many tiny conductors and must therefore be handled with care.

DO NOT STORE THE CAMERA SYSTEM in direct sunlight or at high temperatures.

ENSURE THAT THE CABLE IS NOT COILED TOO TIGHTLY OR FOLDED UP WHILE THE CAMERA IS IN STORAGE. Always arrange the cable in loose coils with a diameter of 6 inches (15 cm) or more.

DO NOT PULL FORCEFULLY ON THE CABLE WHEN UNCOILING. Undesirable 'kinks' may result. Always allow the cable to uncoil gently and loosely.

HANDLE THE SYSTEM WITH CARE, avoiding any strong impacts. When transporting, use the original packaging to prevent damage to the device.

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- This reprocessing step is part of the validated process.
- This reprocessing step may be carried out additionally in order to aid cleaning. The sterilization method is approved in respect of material compatibility, however must be validated by the operator on site in respect of effectiveness.
- x On the grounds of potential damage, this reprocessing step must not be performed.

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7.2.2 General warnings

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

WARNING: Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

CAUTION: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: Damage to product. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be requested from hygiene@karlstorz.com.

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**Reprocessing
(not validated for
US customers)**

CAUTION: National laws and regulations must be observed.

NOTE: The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (art. no. 96216003) can be downloaded or requested by visiting www.karlstorz.com.

7.2.3 Accessories

Accessories required for carrying out reprocessing:

- Brush: 27652

7.2.4 Manual wipe-down disinfection

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixing effect and material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time, wipe the surface with a dry low lint cloth.

7.2.5 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, pre-clean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual pre-cleaning under cold running water.

7.2.6 Manual cleaning

The medical device must be completely immersed in a cleaning solution. To ensure bubble-free wetting, the lumina in particular must be filled. After the necessary exposure time, clean the instrument mechanically with the aid of brushes or sponges. Finally, it must be rinsed with cold water to ensure neutralization.

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7.2.7 Manual disinfection

The medical device must be completely immersed in a disinfectant solution. To ensure bubble-free wetting, the lumina in particular must be filled. At the end of the necessary exposure time, the medical device must be rinsed several times in order to remove all chemical residues. To this end, water of the best possible quality must be used, taking into account the national regulations. Finally, all of the surfaces, joints, openings, channels and lumina are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

7.2.8 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Afterwards, a functional check must be carried out.

NOTE: During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

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7.2.10 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003). The method must be selected taking into account the respective applicable national requirements and in consultation with the device and product manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

**Hydrogen peroxide (H₂O₂)
sterilization – ASP STERRAD®**

WARNING: Medical devices must be sterilized in dismantled condition.

WARNING: Remember that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

CAUTION: For detailed information, please consult the user handbook of the respective device.

NOTE: Refer to the 'STERRAD® Sterility Guide' to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

NOTE: Sterilization is not possible on surfaces which have been lubricated and oiled.

The following STERRAD® sterilization methods have been validated and approved by KARL STORZ for this medical device:

- STERRAD® 50, 100S, 200 without Boosters
- STERRAD® NX® standard cycle
- STERRAD® 100NX® standard cycle

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**Hydrogen peroxide (H₂O₂) sterilization
– STERIS® AMSCO® V-PRO® 1**

Detailed information on the selection of the appropriate cycle for the various device generations is available from the STERIS® manufacturer.

i NOTE: Sterilization is not possible on surfaces which have been lubricated and oiled.

7.2.11 Reprocessing limits

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

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8 Maintenance

8.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be requested from your local representative or from the manufacturer.

8.2 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

CAUTION: Do not attempt to open the camera head or lens. Opening these sealed assemblies affects the camera's soakability and voids all warranties. Do not remove the covers on the camera processor. Warranty is voided unless removed by authorized personnel.

8.3 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this unit.

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8.4 Repair program

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In the case of repairs, please contact your local KARL STORZ subsidiary or authorized dealer.

In Germany you can refer repairs direct to:

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Take-off-Gewerbepark 83
78579 Neuhausen

8.5 Infection prevention

To prevent infections, contaminated devices and instruments must not be returned.

Decontaminate instruments and devices before returning. KARL STORZ returns contaminated instruments/devices to the sender.

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8.6 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet applicable national and local electrical codes, and
- the unit has been used in accordance with its operating instructions at all times.

8.7 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see 'Subsidiaries' section), even during the warranty period.

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**Kompatible
Kamera-Kontrolleinheit**

**Compatible
camera control unit (CCU)**

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**9 Compatible
camera control unit (CCU)**

Item	Cat. No.
TELE PACK +	TP101
C-MAC® Monitor	8403ZX
C-HUB® II	20290320

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**10 Spare parts,
recommended accessories**

Please check your local region for recommended accessory availability.

Article	Order no.
Adaptor, autoclavable, allows telescopes to be changed under sterile conditions	533TVA
Dust cap for camera heads	6349190
Sterile drapes* for KARL STORZ cameras	
Instruction Manual	96286003DER
Instruction manual TELE PACK +	96286007DER
Instruction manual C-MAC® Monitor	96076008DER
Instruction manual C-HUB® II	96206529DER

* For safe and comfortable working with KARL STORZ cameras, we recommend the use of sterile drapes.

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Information on electromagnetic compatibility (EMC)

11 Information on electromagnetic compatibility (EMC)

The H1 camera head TH130 in connection with the camera control unit used is suitable for use within close proximity to an active HF electro-surgical device in professional healthcare facility environment. Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the H1 camera head, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

NOTE: The emission characteristics of this equipment make it suitable for use in industrial areas as well as in hospitals (CISPR 11 Class A) and other professional healthcare environment. If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio communication service. The user might need to take mitigation measures, such as relocating or re-orienting equipment.

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Table 1 – Compliance level for Immunity tests
Guidance and manufacturer's declaration – electromagnetic immunity

The H1 camera head is intended for use in the electromagnetic environment specified below. The user of the H1 camera head should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV Contact discharge ± 15 kV Air discharge	± 8 kV Contact discharge ± 15 kV Air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1 kV ± 1 kV for input and output lines 100 kHz repetition	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage Dip:</u> Dip to 0 % for 1 cycle @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	<u>Voltage Dip:</u> Dip to 0 % for 1 cycle @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the H1 camera head require continued operation during power mains interruptions, it is recommended that the H1 camera head be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50 Hz/60 Hz	30 A/m at 50 Hz/60 Hz	If image distortion occurs, it may be necessary to position the H1 camera head further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3 V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3 V _{rms} on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V _{rms} in ISM bands	3 V _{rms} on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V _{rms} in ISM bands	

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Table 2

Test levels for Proximity fields from RF wireless communications equipment

Test Frequency MHz	Band MHz	Service	Modulation	Immunity Test Level V/m	Compliance level V/m
385	380 – 390	TETRA 400	Pulse modulation 18 Hz	27	27
450	430 – 470	GMRS 460 FRS 460	FM $\pm$ 5 kHz deviation 1 kHz Sine Wave	28	28
710	704 – 787	LTE Band 13 & 17	Pulse modulation 217 Hz	9	9
745					
780					
810	800 – 960	GSM 800/900 TETRA 800 iDEN 820 CDMA 850 LTE Band 5	Pulse modulation 18 Hz	28	28
870					
930					
1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	28	28
1845					
1970					
2450	2400 – 2570	Bluetooth WLAN 802.11 b/g/n RFID 2450 LTE Band 7	Pulse modulation 217 Hz	28	28
5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9	9
5500					
5785					

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Table 3 – Test levels for Radiated and conducted Immunity Tests

Guidance and manufacturer's declaration – electromagnetic immunity

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Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{ms} 150 kHz to 80 MHz	3 V _{ms}	Portable and mobile RF communications equipment should be used no closer to any part of the H1 camera head, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b .
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz Interference may occur in the vicinity of equipment marked with the following symbol: 

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects, and persons.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the H1 camera head is used exceeds the applicable RF compliance level above, the H1 camera head should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the H1 camera head.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 – Emission class and group		
Guidelines and manufacturer's declaration – electromagnetic emissions		
The H1 camera head is intended for use in the electromagnetic environment specified below. The customer or user of the H1 camera head should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The H1 camera head uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class A	The H1 camera head is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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Table 5
Recommended separation distances between portable and mobile RF communications equipment and the H1 camera head

The H1 camera head is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the H1 camera head can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the H1 camera head as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.

HF surgical instrument compatibility

This equipment has been evaluated for compatibility with high-frequency surgical equipment. It has been tested according to IEC 60601-2-2 Annex BB.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@kseas.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kseas.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlántico, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzcube@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Av. Ejército Nacional No. 453 Piso 2,
Colonia Granada, Alameda Miguel Hidalgo
C.P. 11520 Ciudad de México
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Stamvæien 1
1483 Hagan, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
127 39 Skärholmens, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se
KARL STORZ Endoscopy Suomi Oy
Tavastie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestučio st. 59 / Lenktoji st. 27
08124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 67 000
E-Mail: info-it-lv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlystvej 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TU, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopia Nederland B.V.
Displayweg 2
3821 BT Amerstort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N.V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopia France S.A.S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.com

KARL STORZ Endoscopia Ibérica S.A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Croatia d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Cesta v Gorice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Polska Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Interchange 21 Tower, Level 33
399 Sukhumvit Road,
North Khlongtoey, Wattana,
10110 Bangkok, Thailand
Phone: +66 2 660 3669
E-Mail: info-th@karlstorz.com

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, MP Plaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039
E-Mail: info-vietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd.
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch
Room 701A Building 5 & Room 501 Building 7,
No. 3000 Longdong Avenue, Pilot Free Trade Zone,
201203, Shanghai, People's Republic of China
Phone: +86 21 60339888, Fax: +86 21 60339808
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
Room 803-805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
610044, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Sherhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch
Room 02B & 03 & 04A, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 87321281, Fax: +86 20 87321286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: info-asia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: info-asia@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Indonesia
Sinarmas MSG Tower Level 37
Jl. Jend. Sudirman No. Kav. 21
Jakarta Selatan
DKI Jakarta 12920
E-Mail: info-indonesia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info-jp@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
68 Waterloo Road, Macquarie Park NSW 2113
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0895
E-Mail: karlstorz@karlstorz.com.au
www.karlstorz.com

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The subject KARL STORZ products do not incorporate a medicinal substance, nor is it present in the manufacturing process or present in the final device, therefore this section does not apply.

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Recommended transporatation and storage ranges are presented below:

Characteristic	Value
Temperature	From -10 °C to + 60 °C
Relative humidity (non-condensing)	5 %...85 %
Air pressure	500 hPa – 1080 hPa

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This medical device does not harm people and the environment when used for its intended purpose, transportation and storage, subject to compliance with the requirements of regulatory documentation, technical and operational documentation of the manufacturer.

Unused products whose packaging is damaged are class A medical waste of SanPiN 2.1.3684-21. The products used are medical waste of class B in accordance with SanPiN 2.1.3684-21.

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Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications or repairs have been performed by a person or persons duly authorized by KARL STORZ,
- all electrical installations at the location in which the device is connected and operated meet the applicable laws and standards, and
- the device has been used in accordance with its operating instructions at all times.

Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local, even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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

videoscope can be used in combination with the C-HUB II and an external monitor.

2. C-MAC Monitor 8403ZX


Registration certificate - RZN 2019/9503 31.12.2019

C-MAC Monitor for CMOS Endoscopes, screen size 7" with 1280 x 800 pixel resolution, two camera inputs, a USB and a HDMI port, optimized user interface, video and image capture in real time on SD card, playback of recorded video clips and still images, data transfer from SD card to USB flash drive possible, splash-proof according to IP54, suitable for wipe disinfection, shockresistant ABS plastic housing, intelligent power management with rechargeable Li-Ion batteries, without mains adaptor and accessories.

3. Compact endoscopic video system TELE PACK +

The product is in the process of registration.

The TELE PACK + is a device which combines the following functions in a compact system:

- Illumination
- Image display
- Image processing and
- Extensive treatment documentation.

Thanks to its novel design, which includes a complete glass front for optimal cleaning and as a touchscreen, it combines state-of-the-art technology with the unique image quality of KARL STORZ camera heads and videoscopes.

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4. Video camera control unit TELECAM C3 TC100RU

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Telecam C3 is one box with two different internal modules Connect Module and Front End Module. The Connect Module works in conjunction with the Front End Module to display a video image on an external medical monitor. Telecam C3 has receptacles to which a camera head is attached to collect the source image, then relay it to the Connect Module.

Compatible endoscopes:

The coupling mechanism allows all KARL STORZ rigid and flexible standard endoscopes to be connected.

These control units are at the same time a power supply for the camera heads and together with the camera heads, light source and endoscopes are considered the ME System. Only the camera control units specified in this Addendum should be used with the system.

The rated voltage of the supply mains of the camera control units compatible with the product, which are the power source of the camera heads, does not exceed 250 V.

Input electrical characteristics of control units

Characteristic	Power supply	Operating frequency	Power input
TELE PACK+	100-240 VAC	50/60 Hz	140 VA
C-MAC	12 VDC	-	18 VA
C-HUB II	5 VDC		3.8 VA
TELECAM C3	100-240 VAC	50/60 Hz	100-135 VA

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• Lowering the status of contamination to unavoidable residual risks

The KARL STORZ material compatibility list is your guide for the reprocessing of rigid and flexible endoscopes as well as endoscopic instruments and provides the best possible value preservation of your investments. The chemicals listed here have been tested and approved by KARL STORZ with regard to material compatibility.

A differentiation is made between manual cleaning and disinfection/ultrasonic cleaning and machine cleaning and disinfection for which various chemical agents are used. Material compatibility approval is based on tests on various processing chemicals. Therefore, any possible chemical cross-reactions due to insufficient neutralization/rinsing between process steps or incorrect dosage/exposure time can not be ruled out. This particularly applies to subsequent low-temperature sterilization as this process also uses chemical agents.

To ensure efficacious reprocessing conforming to legal and statutory requirements, please observe user manuals for each specific product and the detailed user manual “Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments”.

Note:

Current situation

KARL STORZ SE & Co. KG does currently not test any reprocessing chemicals for material compatibility. This program is stopped and will be replaced by a new concept which is actually under development. At this point, no further chemicals will be added to the list of approved detergents any more until the new program is launched.

In the meantime, the options for detergent manufacturers are to either conduct individual tests on flexible endoscopes and standard instruments to confirm material compatibility, or to compare their chemicals with one that is already listed as approved for material compatibility and base a confirmation on this assessment where applicable.

1. Manual cleaning and disinfection/ultrasonic cleaning

The substances listed here have been approved and tested by KARL STORZ with regard to material compatibility. This list is not a statement concerning the microbiological effectiveness of the products. Please contact the chemical manufacturer regarding microbiological effectiveness.

The following KARL STORZ instruments **are not** suitable for **ultrasonic cleaning**:

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- Rigid telescopes
- Flexible and semiflexible fiberscopes
- Other instruments with optical components made of glass, with chips or with electronic components

1.1 General warnings

Warning:

During all cleaning and disinfection procedures on contaminated instruments, please observe the occupational safety guidelines of the employers' liability insurance association (Berufsgenossenschaft) or equivalent organization.

Warning:

Gigasept FF from the chemical manufacturer Schulke & Mayr GmbH at a concentration of 8% and a contact time of 60 minutes is not approved for fiberscopes.

Warning:

In the case of flexible endoscopes, first vent and then remove the leakage tester so that the pressure compensation valve closes again.

Warning:

In the case of flexible endoscopes, the pressure compensation cap must be removed from the vent connection during liquid immersion/chemical disinfection.

Caution:

Please follow the instructions "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments." These provide a detailed description of cleaning, disinfection, and sterilization procedures.

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Caution:

Damage to materials: Cleaning the outer sheath of flexible fiberscopes with ammonium chlorides will damage the material, particularly at the hose nozzle.

Caution:

Cleaning all of the materials with alcohol-based cleaning agents can lead to damage to the products as well as contamination becoming fixed to the surface of the products.

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! Caution:

Please note that when preparing and using the solutions, the **chemical manufacturers' instructions on concentration and contact time must be strictly followed.** Excessively long immersion or incorrect concentrations may cause damage. Please contact the chemical manufacturer regarding microbiological effectiveness.

! Caution:

Using SODIUM HYDROXIDE SOLUTION may lead to surface alterations in aluminum parts (even if coated), plastics, and soldered joints, and it may reduce the service life of the instruments.

! Caution:

Before placing the flexible endoscope in the solution, check leak tightness since any liquid penetration may cause damage.

! Caution:

Regardless of the selected cleaning and disinfection procedure, **flexible endoscopes** should **not** be heated to temperatures **above 60 °C.**

! Caution:

Endoscopes should not be cleaned in an ultrasound bath.

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i Note:

National laws and regulations must be observed.

In the case of semi-critical medical devices that are used in a non-sterile condition, disinfectants with full virucidal efficacy should be used.

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① Note:

The use of cleaning and disinfectant agents with pH values outside the neutral range can cause color-changes in color-anodized components. These changes are cosmetic and do not affect the functionality of the components.

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Bandelin electronic GmbH	STAMMOPUR GR	✓	-	✓	✓	✓	✓	-	-
Bandelin electronic GmbH	STAMMOPUR R	✓	-	✓	✓	✓	✓	-	-
Bochemie a.s.	Chirosan	-	✓	-	✓	✓	✓	-	-
Bochemie a.s.	Chirosan plus	✓	✓	-	✓	✓	✓	-	-
Bochemie a.s.	Discleen Endo AF	-	✓	-	✓	✓	✓	-	-
Bochemie a.s.	Discleen Endo PAA	✓	✓	-	✓	✓	✓	✓	✓
Bochemie a.s.	Discleen enzyme	✓	-	✓	✓	✓	✓	✓	✓
Bochemie a.s.	Discleen Extra	✓	✓	-	✓	✓	✓	✓	✓
Bode Chemie GmbH	Bodedex® forte	✓	-	✓	✓	✓	✓	✓	✓
Bode Chemie GmbH	Bomix® plus	✓	✓	✓	✓	✓	✓	✓	✓
Bode Chemie GmbH	Dismoclean® twin basic/ Dismoclean® twin zyme	✓	-	✓	✓	✓	-	-	-
Bode Chemie GmbH	Korsolex® basic	-	✓	-	✓	✓	✓	✓	✓
Bode Chemie GmbH	Korsolex® extra	-	✓	-	✓	✓	✓	✓	✓
Bode Chemie GmbH	Korsolex® FF	-	✓	-	-	-	-	✓	✓
Bode Chemie GmbH	Korsolex® med AF	-	✓	-	✓	✓	✓	-	-
Bode Chemie GmbH	Korsolex® PAA	-	✓	-	✓	✓	✓	✓	✓
Borer Chemie AG	deconex® 50 FF	-	✓	-	✓	✓	-	-	-

Borer Chemie AG	deconex® 53 INSTRUMENT	✓	✓	✓	✓	✓	✓	-	-
Borer Chemie AG	deconex® 53 PLUS	✓	✓	✓	✓	✓	✓	✓	✓
Borer Chemie AG	deconex® 54 SPORCIDE	-	✓	-	✓	✓	✓	✓	✓
Borer Chemie AG	deconex® HLD PA/PA20	-	✓	-	✓	✓	✓	✓	✓
Borer Chemie AG	deconex® POWER ZYME	✓	-	✓	✓	✓	✓	✓	✓
Borer Chemie AG	deconex® TWIN PH10 deconex® TWIN ZYME	✓	-	✓	✓	✓	-	-	-
Calgon Vestal Laboratories	Klenzyme	✓	-	✓	-	-	-	✓	✓
DuPont	PeraSafe	-	✓	-	✓	✓	✓	-	-
DuPont	Rely+On™ PeraSafe®	✓	✓	-	✓	✓	✓	✓	✓
Deppe, Laboratorium Dr. rer. nat.	EndoStar	-	✓	-	✓	✓	✓	✓	✓
Deppe, Laboratorium Dr. rer. nat	InstruPlus®	✓	✓	✓	✓	✓	✓	✓	✓
Deppe, Laboratorium Dr. rer. nat	InstruPlus N	-	✓	✓	✓	✓	✓	-	-
Deppe, Laboratorium Dr. rer. nat	InstruStar®	✓	✓	✓	✓	✓	-	-	-
Deppe, Laboratorium Dr. rer. nat	InstruZym®	✓	-	✓	✓	✓	✓	✓	✓
DERCLEAN	DERDEVICE PLUS A2	-	✓	-	✓	✓	✓	-	-
DERCLEAN	DERDEVICE PLUS OPA	-	✓	-	✓	✓	✓	-	-
Detrox	Detro Forte®	-	✓	-	✓	✓	-	-	-

Detrox	Detro OPA®	-	✓	-	✓	✓	-	-	-
Detrox	Detro Plus®	-	✓	-	✓	✓	✓	-	-
Detrox	Detro Plus OPA®	-	✓	-	✓	✓	✓	-	-
Ecolab Deutschland GmbH	Sekusept™ aktiv	✓	✓	-	✓	✓	✓	✓	✓
Ecolab Deutschland GmbH	Sekusept™ Cleaner	✓	-	✓	✓	✓	-	✓	✓
Ecolab Deutschland GmbH	Sekusept™ easy	-	✓	-	✓	✓	✓	✓	✓
Ecolab Deutschland GmbH	Sekusept™ Extra N	✓	✓	✓	✓	✓	✓	-	-
Ecolab Deutschland GmbH	Sekusept™ forte	-	✓	✓	✓	✓	✓	✓	✓
Ecolab Deutschland GmbH	Sekusept™ MultiEnzyme	✓	-	✓	✓	✓	✓	✓	✓
Ecolab Deutschland GmbH	Sekusept™ PLUS	✓	✓	✓	✓	✓	✓	✓	✓
Ecolab Deutschland GmbH	Sekusept™ Pulver CLASSIC	✓	✓	-	✓	✓	✓	✓	✓
Esteer Pharma GmbH	ULTRA- DESMIT® AF	✓	✓	-	✓	✓	✓	-	-
Esteer Pharma GmbH	ULTRASEPTIN® activ	✓	✓	-	✓	✓	✓	-	-
Esteer Pharma GmbH	ULTRASEPTIN® Classic	✓	✓	✓	✓	✓	✓	-	-
Esteer Pharma GmbH	ULTRASEPTIN® PAA	✓	✓	-	✓	✓	✓	✓	✓
Fresenius AG	Afid	✓	✓	-	✓	✓	✓	-	-
Fresenius AG	Afid plus	-	✓	-	✓	✓	✓	-	-

Fresenius AG	Sporizid	-	✓	-	-	-	-	✓	-
Fresenius AG	Sporizid FF	✓	✓	-	✓	✓	✓	✓	✓
Heck Dental GmbH	Ventisept M Plus Neu	✓	✓	-	✓	✓	✓	✓	✓
Helvemed	Instrument Forte	✓	✓	✓	✓	✓	✓	✓	✓
HIMED GmbH	EndoSeptol FF	-	✓	-	✓	✓	✓	-	-
HIMED GmbH	InstruCid AF	✓	✓	✓	✓	✓	-	-	-
Holifa Poska Sp. z.o.o	Polsept Holifa	✓	✓	-	✓	✓	✓	-	-
Holifa Poska Sp. z.o.o	Polsept Plus	✓	✓	-	✓	✓	✓	✓	✓
Indeba Industria E Comercio LTDA	Letahdeido	-	✓	-	✓	✓	✓	-	-
INIBSA Laboratorios	INSTRUNET F.A. CONCEN.	✓	✓	✓	✓	✓	✓	-	-
Jose Collado S.A.	Darodor® 4000 Liquido	✓	-	-	✓	✓	✓	-	-
Jose Collado S.A.	Darodor® 9000	-	✓	-	✓	✓	✓	-	-
Jose Collado S.A.	Darodor® Sinaldehyd 2000	-	✓	-	✓	✓	✓	-	-
Lonza AG	Formulierung ID 50	-	✓	-	✓	✓	✓	-	-
Lysoform Dr. Rosemann GmbH	Almyrol®	✓	✓	✓	✓	✓	✓	-	-
Lysoform Dr. Rosemann GmbH	Desoform®	✓	✓	✓	✓	✓	✓	-	-

Lysoform Dr. Rosemann GmbH	Lysoformin® 2000	✓	✓	✓	✓	✓	✓	-	-
Lysoform Dr. Rosemann GmbH	Lysoformin® 3000	✓	✓	-	✓	✓	✓	-	-
Merz Hygiene GmbH	Edisonite Super	✓	-	✓	✓	✓	✓	-	-
Merz Hygiene GmbH	Mucadont-IS	-	✓	✓	✓	✓	✓	-	-
Merz Hygiene GmbH	Mucocit-T neu	✓	✓	✓	✓	✓	✓	✓	✓
Merz Hygiene GmbH	Mucadont- Zymaktiv	✓	-	✓	✓	✓	✓	✓	✓
METASYS Medizintechnik GmbH	ID 50	-	✓	-	✓	✓	✓	-	-
METASYS Medizintechnik GmbH	Green & Clean ID	✓	✓	✓	✓	✓	-	-	-
Metrex Research Corporation	Coldspor	-	✓	-	-	-	-	✓	✓
Metrex Research Corporation	EmPower™	✓	-	-	-	-	-	✓	✓
Metrex Research Corporation	MetriCide™	-	✓	-	-	-	-	✓	✓
Metrex Research Corporation	MetriCide™ OPA Plus Solution	-	✓	-	✓	✓	✓	✓	✓
Metrex Research Corporation	MetriZyme™	✓	-	-	-	-	-	✓	✓
Miltex, Inc.	EZ – Zyme	✓	-	-	-	-	-	✓	✓
Minntech International	Adaspor®	-	✓	-	✓	✓	✓	-	-
Mir Dezinfektsii	MIROXIDE® – 2000	-	✓	-	✓	✓	✓	-	-
Necmkimya	BIORAD GDA 2%	-	✓	-	✓	✓	-	-	-

Necmkimya	BIORAD Rayd Guanide	✓	✓	✓	✓	✓	-	-	-
Nova-Praxis-Hygiene	NOVA Biguades AF	✓	✓	✓	✓	✓	✓	-	-
PRISMAN GmbH	POW(D)ER Konz. IDP 700	✓	✓	-	✓	✓	✓	✓	✓
PRISMAN GmbH	IDP-ic 70 POWDER	✓	✓	-	✓	✓	✓	✓	✓
Profarma	Ajatin® Extra Strong	-	✓	-	✓	✓	✓	✓	✓
Promagent AB	Wavacide	-	✓	-	✓	✓	✓	✓	✓
Rheosol Wachendorff – Chemie GmbH	RHEOSEPT-ID aktiv Granulat	✓	✓	✓	✓	✓	✓	✓	✓
Rheosol Wachendorff – Chemie GmbH	RHEOSEPT-ID plus+	✓	✓	✓	✓	✓	✓	✓	✓
Ruhof Corporation	Endozime®	✓	-	-	-	-	-	✓	✓
Schumacher, Dr. GmbH	DESCOTON extra	✓	✓	-	✓	✓	✓	✓	✓
Schumacher, Dr. GmbH	DESCOTON forte	-	✓	-	✓	✓	✓	✓	✓
Schumacher, Dr. GmbH	DESCOTON 2% GDA	-	✓	-	✓	✓	✓	✓	✓
Schumacher, Dr. GmbH	INSTRU PLUS flussig	✓	-	✓	-	-	-	✓	✓
Schumacher, Dr. GmbH	MULTIZYME FORTE	✓	-	✓	✓	✓	✓	✓	✓
Schumacher, Dr. GmbH	PERFEKTAN® ENDO	✓	✓	✓	✓	✓	✓	✓	✓
Schumacher, Dr. GmbH	PERFEKTAN® NEU	✓	✓	✓	✓	✓	✓	✓	✓
Schumacher, Dr. GmbH	PERFEKTAN® TB	✓	✓	✓	✓	✓	-	-	-

Schumacher, Dr. GmbH	PLURAZYME	✓	-	✓	✓	✓	✓	✓	✓
Schulke & Mayr GmbH	gigasept AF® forte	✓	✓	-	✓	✓	✓	✓	✓
Schulke & Mayr GmbH	gigasept FF® (neu)	-	✓	-	✓	✓	✓	-	-
Schulke & Mayr GmbH	gigasept® instru AF	✓	✓	✓	✓	✓	-	-	-
Schulke & Mayr GmbH	gigasept® med forte	✓	✓	✓	✓	✓	✓	✓	✓
Schulke & Mayr GmbH	gigasept® PAA concentrate	-	✓	-	✓	✓	✓	✓	✓
Schulke & Mayr GmbH	gigazyme®	✓	-	✓	✓	✓	✓	✓	✓
Schulke & Mayr GmbH	Lysetol® V	-	✓	-	✓	✓	✓	-	-
Schulke & Mayr GmbH	mikrozyd® PAA	-	✓	-	✓	✓	✓	✓	✓
Sirmaxo Chemicals Pvt. Ltd.	Sanidex OPA	-	✓	-	✓	✓	✓	-	-
Sirmaxo Chemicals Pvt. Ltd.	Saniscope	-	✓	-	✓	✓	✓	✓	✓
Sirmaxo Chemicals Pvt. Ltd.	Sanizyme	✓	-	-	✓	✓	✓	✓	✓
Steril-4 s.r.l.	STERIL-C	-	✓	-	✓	✓	✓	-	-
STERIS Corporation	EnzyCare 2	✓	-	✓	✓	✓	✓	-	-
STERIS Corporation	Revital-Ox® Resert™ HLD	-	✓	-	✓	✓	✓	✓	✓
STERIS Corporation	Revital-Ox™ Resert™ HLD Chemosterilant	-	✓	-	✓	✓	✓	✓	✓
STERIS Corporation	Resert™ XL HLD	-	✓	-	✓	✓	✓	✓	✓

Tristel Solutions Ltd.	Medistel	-	✓	✓	✓	✓	✓	-	-
Tristel Solutions Ltd.	Enzystel	✓	-	✓	✓	✓	✓	-	-
Tristel Solutions Ltd.	Stella™ System mit: Tristel™ Clean for Stella Tristel Fuse™ for Stella	✓	✓	-	-	-	✓	✓	-
Tristel Solutions Ltd.	Perastel 3 Day	-	✓	-	✓	✓	✓	✓	✓
Tristel Solutions Ltd.	Tristel Wipes System	✓	✓	-	✓	-	✓	✓	✓
Wavi Energie System INC.	Wavicide-01	-	✓	-	-	-	-	✓	✓
Weigert, Dr. GmbH & Co.	neodisher endo® CLEAN	✓	-	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher endo® DIS active	✓	✓	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® LM 2	✓	-	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® MediClean	✓	-	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® MediClean forte	✓	-	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® MediZym	✓	-	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® PreStop	✓	-	-	✓	✓	-	-	-
Weigert, Dr. GmbH & Co.	neodisher® Septo 3000	-	✓	✓	-	-	-	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® Septo Active	-	✓	✓	✓	✓	-	-	-
Weigert, Dr. GmbH & Co.	neodisher® SeptoClean*	✓	✓	✓	✓	✓	-	-	-

Weigert, Dr. GmbH & Co.	neodisher® Septo MED	✓	✓	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® Septo PreClean	✓	✓	✓	✓	✓	✓	✓	✓
Whiteley Industries Pty. Ltd.	Aidal® Plus	-	✓	-	✓	✓	✓	-	-
Whiteley Industries Pty. Ltd.	Matrix®	✓	-	-	✓	✓	✓	✓	✓

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! Caution:

To minimize stress on materials, avoid repeated switching between processing procedures.

! Caution:

Flexible endoscopes should never be steam sterilized (autoclaved).

! Caution:

The pressure compensation cap must be removed from the vent connection during liquid immersion/chemical disinfection.

! Caution:

Please note that when preparing and using the solutions, the **chemical manufacturers' instructions on concentration and contact time must be strictly followed**. Excessively long immersion or incorrect concentrations may cause damage. Please contact the chemical manufacturer regarding microbiological effectiveness.

! Caution:

Regardless of the selected cleaning and disinfection procedure, **flexible endoscopes** should **not** be heated to temperatures **above 60 °C**.

! Caution:

The use of rinse aids is not approved for plastics. Rinse aids have a negative effect on plasticizers, which may lead to breakage.

! Caution:

The product thermosept RKF forte from Schülke & Mayr GmbH is not suitable for aluminum materials.

i Note:

The C-MAC® video laryngoscope series (840XXX) is considered to be a video endoscope.

i Note:

For detailed information on processing, please refer to the respective instructions for use.

i Note:

Please contact the chemical manufacturer regarding the range of action of the chemical agents used.

i Note:

The use of cleaning and disinfectant agents with pH

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values outside the neutral range can cause color-changes in color-anodized components. These changes are cosmetic and do not affect the functionality of the components.

Note:

When using cleaning and disinfecting agents, chemicals should be carefully coordinated in order to avoid any chemical interactions leading to material damage. If you have any questions or doubts, please contact the chemical manufacturer.

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Bode Chemie GmbH	Dismoclean® 28 Alka one	✓	-	✓	✓	✓	✓	-	-
Bode Chemie GmbH	Dismoclean® twin basic/ Dismoclean® twin zyme	✓	-	✓	✓	✓	-	-	-
Bode Chemie GmbH	Korsolex® Endo- Cleaner	✓	-	✓	✓	✓	✓	✓	✓
Bode Chemie GmbH	Korsolex® Endo- Disinfectant	-	✓	✓	✓	✓	✓	✓	✓
Borer Chemie AG	deconex® ALKAMILD	✓	-	-	✓	✓	✓	-	-
Borer Chemie AG	deconex® ENDOMATIC	-	✓	-	-	-	✓	✓	-
Borer Chemie AG	deconex® POWER ZYME	✓	-	✓	✓	✓	✓	✓	✓
Borer Chemie AG	deconex® TWIN PH10 deconex® TWIN ZYME	✓	-	✓	✓	✓	-	-	-
B. Braun Medical AG	Helimatic® Cleaner alkaline	✓	-	-	✓	✓	✓	-	-
B. Braun Medical AG	Helimatic® Cleaner enzymatic	✓	-	-	✓	✓	✓	✓	✓
B. Braun Medical AG	Helimatic® Cleaner neutral	✓	-	-	✓	✓	-	-	-
B. Braun Medical AG	Helimatic® Disinfectant	-	✓	-	✓	✓	✓	✓	✓
Deppe, Laboratorium Dr. rer. nat.	Endomat Plus	-	✓	-	✓	✓	✓	✓	✓
Deppe, Laboratorium Dr. rer. nat.	EndomatPlus viruguard	-	✓	-	✓	✓	✓	✓	✓
Deppe, Laboratorium Dr. rer. nat.	Labomat E	✓	-	-	✓	✓	✓	✓	✓
DiverseyLever	Sumatox E	✓	-	✓	✓	✓	✓	-	-

Ecolab Deutschland GmbH	MetalClean Plus	✓	✓	-	✓	✓	✓	-	-
Ecolab Deutschland GmbH	Sekumatic™ FD	-	✓	-	✓	✓	✓	-	-
Ecolab Deutschland GmbH	Sekumatic™ FR	✓	-	-	✓	✓	✓	-	-
Ecolab Deutschland GmbH	Sekumatic™ FRE	✓	-	-	✓	✓	✓	-	-
Ecolab Deutschland GmbH	Sekumatic™ MultiClean	✓	-	-	✓	✓	✓	-	-
IMS	Adapspor® Ready to Use	-	✓	✓	✓	✓	✓	-	-
Medisafe UK, Ltd.	3E-Zyme/HS-Zyme	✓	-	✓	✓	✓	✓	-	-
MEDIVATORS Inc.	Rapicide®	-	✓	-	✓	✓	✓	✓	✓
MEDIVATORS Inc.	Rapicide® PA	-	✓	-	✓	✓	✓	✓	✓
Merz Consumer Care GmbH	mucapur® AF	✓	-	✓	✓	✓	✓	-	-
Merz Consumer Care GmbH	mucapur® ED	-	✓	✓	✓	✓	✓	✓	✓
Merz Consumer Care GmbH	mucapur® ED plus	-	✓	✓	✓	✓	✓	✓	✓
Merz Consumer Care GmbH	mucapur® ER	✓	-	✓	✓	✓	✓	✓	✓
Merz Consumer Care GmbH	mucapur® ER plus	✓	-	✓	✓	✓	✓	✓	✓
Ruhof Corporation	Endozime® AW	✓	-	✓	✓	✓	✓	-	-
Schumacher, Dr. GmbH	THERMOTON® CLEANER	✓	-	-	✓	✓	✓	-	-
Schumacher, Dr. GmbH	THERMOTON® ENDO	-	✓	✓	✓	✓	✓	✓	✓

Schumacher, Dr. GmbH	THERMOTON [®] DESIN- FEKTANT	-	✓	-	✓	✓	✓	✓	✓
Schumacher, Dr. GmbH	THERMOTON [®] NR	✓	-	-	✓	✓	✓	✓	✓
Schulke & Mayr GmbH	thermosept [®] alka clean	✓	-	-	✓	✓	✓	-	-
Schulke & Mayr GmbH	thermosept [®] alka clean forte	✓	-	-	✓	✓	✓	-	-
Schulke & Mayr GmbH	thermosept [®] DK	-	✓	-	✓	✓	✓	-	-
Schulke & Mayr GmbH	thermosept [®] ED	-	✓	-	✓	✓	✓	✓	✓
Schulke & Mayr GmbH	thermosept [®] ER	✓	-	-	✓	✓	✓	✓	✓
Schulke & Mayr GmbH	thermosept [®] PAA base thermosept [®] PAA additive	-	✓	-	✓	✓	✓	✓	✓
Schulke & Mayr GmbH	thermosept [®] RKF forte	✓	-	-	✓	✓	✓	-	-
Schulke & Mayr GmbH	thermosept [®] RKN- zym	✓	-	-	✓	✓	✓	-	-
Schulke & Mayr GmbH	thermosept [®] X-tra	✓	-	-	✓	✓	-	-	-
STERIS Corporation	Instru-Klenz	✓	-	✓	✓	✓	✓	-	-
Weigert, Dr. GmbH & Co.	neodisher endo [®] CLEAN	✓	-	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher endo [®] SEPT GA	-	✓	-	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher endo [®] SEPT PAC	-	✓	-	-	-	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher [®] DuoClean	✓	-	-	✓	✓	-	-	-

Weigert, Dr. GmbH & Co.	neodisher® FA	✓	-	-	✓	✓	-	-	-
Weigert, Dr. GmbH & Co.	neodisher® MediClean	✓	-	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® MediClean forte	✓	-	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® MediZym	✓	-	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® OXIVARIO	✓	-	-	✓	✓	-	-	-
Weigert, Dr. GmbH & Co.	neodisher® SeptoClean*	✓	✓	✓	✓	✓	-	-	-
Weigert, Dr. GmbH & Co.	neodisher® Septo DN	-	✓	-	✓	✓	✓	✓	✓

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✓ Suitable chemical agents

- Chemicals not tested

Cleaning agents that are marked with*:

These cleaning agents are proven to be effective against prions

① Note:

If agents other than those listed above have been used, any arising damage is not covered by KARL STORZ.

Here with KARL STORZ SE & Co. KG confirms the material compatibility of the following chemicals for processing endoscopes and endoscopic instruments.

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