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- Incorrect application of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

- Abdominal puncture sites 7 mm or greater (for the introduction of auxiliary instrumentation) may be a source of herniation.

- During monopolar electrosurgery, inadvertent burns can occur when the appropriate patient return path is obstructed or as a result of capacitive coupling.

Risk of injury and damage to the products:

- Failure to observe and follow video endoscope instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals and always follow the instructions given precisely. Check the functioning of the products used in combination.

- The use of application parts beyond one's field of vision bears a risk of unintended tissue injury and damage to accessories. During use always hold the application parts of the active electrode as well as of the laser and other energy conducting equipment within your field of vision.

- The high light energy could heat up the distal end, the light connection, as well as adjacent components and material in front of the light emission window. This can lead to burns in the patient, the user or the surgical equipment.

- Avoid direct material contact with the distal end and the light connection of the endoscope and the light cable.

- Do not leave the endoscope and the light conductor on the patient or in direct contact with surgical equipment.

- Always select a light output for the cold light source that is as low as possible to barely achieve optimal illumination of the surgical field.

Risk of fatal injury:

In the case of patients with cardiac pacemakers it must be ensured that the proximal end (furthest from the patient) of the TIPCAM 1 Rubina is not positioned less than 2 cm from the pacemaker. The device generates a magnetic field at the handle, which can cause cardiac pacemakers or other electronic devices to malfunction.

Risk of electric shock:

Do not attempt to open the TIPCAM 1 Rubina. Refer servicing to qualified

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Outer diameter	10,3 mm ± 0,1 mm
Instrument channel width	N/A, no channel
Cable length	2,9m±0,1m
Direction of view	0° ± 3° 30° ± 3°
Field of view (in 2D WL and 3D WL)	82° ± 5°
Field of view (in 2D Auto Leveling)	71,3° ± 3°
Depth of focus	30 ... 200 mm
Sensor type	CMOS
Frame rate	50/60 Гц
Resolution	3840 x 2160
Weight	420 g ± 30 g (excluding the video cable)

Key accessories, medical devices or non-medical devices, to be used in combination with medical device submitted for registration

List of compatible devices

№ RC	Name
CCU	
PY P3H 2015/3274	Endoscopic camera control unit IMAGE1 S CONNECT II
PY P3H 2015/3274	Endoscopic camera control unit IMAGE1 S 4U-LINK
PY P3H 2015/3274	Endoscopic camera control unit IMAGE1 S CONNECT
PY P3H 2015/3274	Endoscopic camera control unit IMAGE1 S H3-LINK
Light guides	
P3H 2016/4855	Fiber optic light-guide cable, Ø 4.8 mm, length 300 cm, extra heat resistant (495TIP)
P3H 2016/4855	Fiber optic light guide cable, Ø 4.8 mm, length 250 cm, extra heat resistant, safety interlock (495 NCSC)
Sources of light	
P3H 2021/15229	Endoscopic cold light source Power LED Rubina

Information about the medicines included in the medical device

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Dispose of the system in accordance with local regulations for the disposal of such products.

If you want to stop using the product, please note that state municipal regulations may contain special rules for the disposal of this equipment. Contact Karl Storz Technical Support to ensure compliance with these legal provisions and to avoid environmental damage that could result from the disposal of your system.

The packaging material must be disposed of in accordance with the environmental requirements of national legislation.

The system and accessories must be disposed of at the end of their useful life in accordance with applicable regional, state, and local codes.

For more information, contact your local service representative or Karl Storz regional office.

Products coming into contact with blood and/or biological fluids are class B medical devices under medical wastes classification and may be recycled and / or disposed in accordance with applicable national regulations.

Unused products (not having contact with blood and/or biological fluids) belong to class A medical devices under medical waste classification and may be recycled and/or disposed in accordance with applicable national regulations.

List of international regulatory documents/standards which the medical device complies with

Standard №	Title
ISO 14971:2019	Medical devices — Application of risk management to medical devices
ISO 13485:2016	Medical devices — Quality management systems — Requirements for regulatory purposes
ISO 15223-1:2016	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements

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IEC 62366-1:2015	Medical devices — Part 1: Application of usability engineering to medical devices
ISO 8600-1:2015	Endoscopes — Medical endoscopes and endotherapy devices — Part 1: General requirements
ISO 8600-3:2019	Endoscopes — Medical endoscopes and endotherapy devices — Part 3: Determination of field of view and direction of view of endoscopes with optics
ISO 8600-4:2014	Endoscopes — Medical endoscopes and endotherapy devices — Part 4: Determination of maximum width of insertion portion
ISO 10993-1:2018	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
ISO 10993-5:2009	Biological evaluation of medical devices - Part 5: Tests for <i>in vitro</i> cytotoxicity
ISO 10993-10:2010	Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
ISO 10993-11:2017	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity
ISO 10993-12:2012	Biological evaluation of medical devices - Part 12: Sample preparation and reference materials
ISO 10993-17:2002	Biological evaluation of medical devices - Part 17: Establishment of allowable limits for leachable substances
ISO 10993-18:2020	Biological evaluation of medical devices — Part 18: Chemical characterization of medical device materials within a risk management process
ISO 17664:2017	Processing of health care products — Information to be provided by the medical device manufacturer for the processing of medical devices
ISO 17665-1:2006	Sterilization of health care products — Moist heat — Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices
ISO 14937:2009	Sterilization of health care products — General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices
ISO 15883-1:2006	Washer-disinfectors — Part 1: General requirements, terms and definitions and tests
IEC 60601-1:2005	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance

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contact the manufacturer and / or authorized representative of the manufacturer.

Manufacturer:

Karl Storz SE & Co.KG, Germany
(“Karl Storz CE and Co. KG”, Germany)
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen, Germany
Tel. +49 7461 708-0, +49 7461 708-105
info@karlstorz.com

Authorized representative of the manufacturer:

Limited Liability Company “Karl Storz - Endoscopes Vostok” (LLC “KSHTÉ VOSTOK”)
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: Info-ru@karlstorz.com

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Thank you for placing your trust in the name KARL STORZ. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high quality piece of equipment from KARL STORZ.

This instruction manual is intended to aid you in the proper installation, connection, and operation of this device. All required details and all actions on your part are clearly presented and explained. Please read this manual carefully before proceeding with the equipment. Keep this instruction manual available for ready reference in a convenient location near the equipment.

i NOTE: This IFU is applicable for IMAGE1 S CONNECT® TC200 and for IMAGE1 S CONNECT® II TC201. Only when a special feature of IMAGE1 S CONNECT® TC200 or IMAGE1 S CONNECT® II TC201 is described the model number TC200 or TC201 is mentioned in the IFU to distinguish between CONNECT and CONNECT II. Otherwise the term «IMAGE1 S CONNECT®» is used.

It is recommended to check the suitability of the product for the intended procedure prior to use.

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1 Images of the equipment

1.1 Possible combinations for IMAGE1 S™

The IMAGE1 S™ comprises IMAGE1 S CONNECT® (TC200) or IMAGE1 S CONNECT® II (TC201) and a minimum of one and maximum of three additional LINK modules. The IMAGE1 S™ H3-LINK (TC300), IMAGE1 S™ X-LINK (TC301), IMAGE1 S D3-LINK® module (TC302) and IMAGE1 S™ 4U-LINK module (TC304) are available as LINK modules.

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**2 Controls, displays,
connectors, and their uses**

**2.1 IMAGE1 S CONNECT® TC200
in combination with the
H3-LINK, X-LINK and D3-LINK
modules (see page IV)**

- ① Power switch
- ② White balance
- ③ USB port
- ④ Potential equalization connector
- ⑤ LINK inputs
- ⑥ SCB input/output sockets
(KARL STORZ Communication Bus)
- ⑦ LAN interface
- ⑧ USB interface
- ⑨ HD DVI-D video output sockets
- ⑩ Potential equalization connector
- ⑪ 3G-SDI video output
- ⑫ Power cord socket
- ⑬ LINK output
- ⑭ H3-LINK camera connector
- ⑮ X-LINK camera connector
- ⑯ D3-LINK camera connector

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**2.2 IMAGE1 S CONNECT® II
TC201 in combination with the
4U-LINK, X-LINK and D3-LINK
modules (see page V)**

- ① Power switch
- ② White balance
- ③ USB port
- ④ Potential equalization connector
- ⑤ LINK inputs
- ⑥ SCB input/output socket
(KARL STORZ Communication Bus)
- ⑦ LAN interface
- ⑧ USB interface
- ⑨ HD DVI D video output socket
- ⑩ Potential equalization connector
- ⑪ 12G/3G-SDI video output
- ⑫ Power cord socket
- ⑬ LINK output
- ⑭ 4U-LINK camera connector
- ⑮ X LINK camera connector
- ⑯ D3-LINK camera connector
- ⑰ DisplayPort video output sockets
- ⑱ EXT SYNC connector
(for future use)

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Non-ionizing electromagnetic radiation

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

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

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1 Images of the equipment	IV
1.1 Possible combinations for IMAGE1 S™	IV
2 Controls, displays, connectors, and their uses	VI
2.1 IMAGE1 S CONNECT® TC200 in combination with the H3-LINK, X-LINK and D3-LINK modules (see page IV)	VI
2.2 IMAGE1 S CONNECT® II TC201 in combination with the 4U-LINK, X-LINK and D3-LINK modules (see page V)	VII
3 Symbols employed	VIII
3.1 Symbols for operation	VIII
3.2 Label and packaging symbols	IX
4 General information	4
4.1 Description of the devices	4
4.2 Property rights	4
5 Safety instructions	5
5.1 Explanation of warnings and cautions	5
5.2 General safety information	6
5.3 Indications for use	10
5.3.1 Intended use	10
5.4 Contraindications	10
5.5 Patient profile	11
5.6 User profile (specialist physician)	11
5.7 User profile (assistant)	12
5.8 Safety precautions at the site of installation	13
5.9 Safety precautions when operating the device	13
6 Initial operation	14
6.1 Basic equipment	14
6.1.1 TC200 IMAGE1 S CONNECT®	14
6.1.2 TC201 IMAGE1 S CONNECT® II	14
6.1.3 TC300 IMAGE1 S™ H3-LINK	14
6.1.4 TC301 IMAGE1 S™ X-LINK	14
6.1.5 TC302 IMAGE1 S D3-LINK®	14
6.1.6 TC304 IMAGE1 S™ 4U-LINK	14
6.2 Unpacking the equipment	15
6.3 Preparing the device for use	15
6.3.1 Installation	15
6.3.2 Connecting the potential equalization	15

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6.3.3	Connecting the IMAGE1 S™ LINK module with IMAGE1 S CONNECT®.....	15
6.3.4	Connecting SCB.....	16
6.3.5	Establishing a KS HIVE connection.....	16
6.3.6	Connecting the monitor.....	16
6.3.7	Connecting USB devices.....	17
6.3.8	Connecting the power supply.....	17
6.3.9	Switching on the IMAGE1 S CONNECT®.....	17
6.3.10	Connecting USB storage media.....	18
6.3.11	Connecting the camera head/ video endoscope.....	18
6.3.12	Performing the white balance.....	18
6.3.13	Test for proper functioning.....	19
6.3.14	Installation overview (TC200).....	20
6.3.15	Installation overview (TC201).....	21
7	Operating instructions.....	22
7.1	Basic information on operation.....	22
7.1.1	Keyboard and mouse.....	22
7.1.2	Camera head buttons.....	22
7.2	Dashboard.....	22
7.2.1	Description of functions of the symbols used.....	23
7.3	Live Menu.....	27
7.3.1	Description of functions of the symbols used.....	28
7.4	Setup Menu.....	37
7.4.1	General Settings (IMAGE1 S CONNECT®).....	37
7.4.2	Presets.....	42
7.4.3	LINK-specific settings.....	42
7.5	Patient Data Privacy.....	43
8	Reprocessing.....	46
8.1	General safety information.....	46
8.2	Wipe-down disinfection of devices.....	47
9	Service and repair.....	48
9.1	Maintenance and safety check.....	48
9.1.1	Maintenance.....	48
9.1.2	Safety check.....	48
9.2	Servicing and repair.....	49
9.3	Disposal.....	49
9.4	Repairing devices.....	50
9.5	Infection prevention.....	50
9.6	Limitation of liability.....	51
9.7	Warranty.....	51

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4 General information

4.1 Description of the devices

The IMAGE1 S™ takes a unique and innovative approach to system conception. The modular design principle allows the user to create a system which is tailored to his needs. One IMAGE1 S CONNECT™ or IMAGE1 S CONNECT™ II module is at the heart of the IMAGE1 S™ and acts as a communication interface between the image and documentation devices.

For use, one IMAGE1 S CONNECT™ or IMAGE1 S CONNECT™ II plus at least one LINK module is required which receives and forwards the signal. The Light Source Control and S technologies for image optimization provide users with even more comprehensive support. An intuitive user interface assures simple and quick access as well as rapid understanding of menu navigation.

The system offers numerous state-of-the-art functions. These include:

- IMAGE1 S™ visualization modes to assist the user via contrast enhancement and/or more even brightness perception
- Light Source Control for automatic adjustment of the light source
- Integrated image and video capture for simple documentation

Functions are described in greater detail in § 7.4.

4.2 Property rights

Protected by at least one of the following US Patents: 7.821.530; 8.274.559; 7.520.853; 7.212.227; 8.089.509; 8.199.188. For more information, visit: www.karlstorz.com.

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5 Safety instructions

5.1 Explanation of warnings and cautions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

Please read this manual and follow the instructions carefully. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system. To make these words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician is at risk.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

NOTE: A Note indicates special information about operating the device or clarifies important information.

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5.2 General safety information

 **WARNING:** Risk of injury and damage to the products: Failure to observe and follow this instruction manual and all of the instruction manuals of the products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals carefully and always follow the instructions given precisely. Check the functioning of the products used in combination.

 **WARNING:** Test this device prior to each surgical procedure to ensure that it functions correctly. In the event that the image becomes unusable during surgery, the camera may be disengaged from the endoscope and the procedure continued optically. If this is not possible, it is up to the discretion of the surgeon how best to proceed. Availability of a spare system is recommended.

 **WARNING:** Danger of an electric shock! This device may only be connected to a socket with protective earthing. To reduce the risk of electric shock, keep the device out of reach of patients.

 **WARNING:** Grounding reliability can only be achieved when the equipment is connected to a 'Hospital Only' or 'Hospital Grade' outlet (i.e., approved for use in an operating room environment). The plug and cord should be routinely inspected. Do not use if damaged.

 **WARNING:** Before each use or after a change of viewing modes/settings, the operator should check to ensure the view observed through the endoscope provides a live image (rather than a stored one) and has the correct image orientation.

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 **WARNING:** Consult and follow the storage, transport and operating conditions prior to storing, transporting or operating, respectively. If the device is stored in a warm, humid environment, ensure that it has suitably equilibrated to the operating environment prior to use.

 **WARNING:** To reduce the risk of electrical shock, do not remove cover of unit. Refer servicing to qualified personnel. Removal of covers by unauthorized personnel will void the equipment's warranty.

 **WARNING:** Perceived electrical shocks may be experienced when using the camera system in conjunction with electrosurgical units.

 **WARNING:** High energy radiated light through endoscopes may give rise to high temperatures in front of the light outlet and to the tip of the endoscope. To minimize the risk of burns, always adjust the light source to the minimum illumination intensity necessary to achieve optimum illumination of the endoscopic scene when coupled to the video camera.

 **WARNING:** Refer to the appropriate section of this manual for validated cleaning and disinfection instructions.

 **CAUTION:** Do not store liquids on or above the unit.

 **CAUTION:** In the USA, federal law restricts this device to sale or use by or on the order of a physician.

 **NOTE:** Disposal of these system products, at the end of their useful life, must be in accordance with local regulations.

 **NOTE:** Do not discard as unsorted municipal waste.

 **NOTE:** Discard as electrical/electronic waste; recycle or reuse accordingly.

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 **WARNING:** Under some Single Fault Conditions per the IEC 60601-1 standard, the temperature of the equipment surfaces could get hot and there is a possible risk of a burn if touched.

 **NOTE:** For best thermal performance the TC201/TC200 should be mounted at the top of the equipment stack.

 **WARNING:** No modification of this equipment is allowed.

 **NOTE:** Consult local authorities for reuse/recycle instructions.

 **NOTE:** The presence of maximum mains voltage on SIP/SOPs is not considered a normal condition.

 **NOTE:** Additional equipment connected to medical electrical equipment must be able to prove their compliance with the respective IEC or ISO standards (e.g., IEC 60950 for data processing equipment). Furthermore, all configurations shall comply with the standard requirements for medical systems (see Clause 16 of the 3rd Ed. of IEC 60601-1). Anybody connecting additional equipment to medical electrical equipment is a system configurator and is therefore responsible for the system's compliance with the standard requirements for systems. Please note that local laws take priority over the above-mentioned standard requirements. If in doubt, please consult your local specialist dealer or the technical service department. (Standard / directive references: IEC 60601-1+A1+A2: 6.8.2.c, 19.2.b, 19.2.c, IEC 60601-1: 7.9.2.5, 8.1, 16.2.d, MDD 93/42/EEC: Annex I clause 13.6.c).

 **NOTE:** For use with the product, KARL STORZ recommends the accessories listed in chapter 11. If third-party accessories are used, the user must check the correct operation himself.

 **WARNING:** When performing fluorescence imaging procedures, use suitable endoscopes and recommended light cables.

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WARNING: To ensure safe operation, do not simultaneously touch the device output connectors and the patient.

WARNING: When endoscopes are used with energized endoscopically used accessories, the patient leakage currents may be additive. This is particularly relevant when using type CF devices. In this case, a Type CF endoscopically used accessory should be used in order to minimize the total patient leakage current.

WARNING: During discharge of a defibrillator, the image displayed may take up to 4 seconds to recover and in some cases a power cycle may be required to restore the image.

WARNING: The device should be placed on a flat surface or in a video cart. Do not block air inlets/outlets as overheating may result in fire hazard or device failure during procedure.

WARNING: Failure to perform proper white balance of system may result in color distortion.

WARNING: The operator is responsible to limit physical access to the IMAGE1 S™ from unauthorized persons.

WARNING: In the case of loss of surgical image for greater than 3 seconds, a power cycle may be required to restore the image.

WARNING: When using high frequency surgical equipment, keep the working part of the active electrode in the field of view of the operator to avoid accidental high frequency burns. Avoid contact between the active electrode and metal parts of the endoscope.

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WARNING: Do not connect the USB port of the IMAGE1 S to network devices such as ethernet, Bluetooth, or WIFI dongles. The IMAGE1 S is not intended for such use, and connection of such devices to the USB port can affect its operation.

5.3 Indications for use

5.3.1 Intended use

IMAGE1 S™ is a camera control unit (CCU) for use with camera heads or video endoscopes for the visualization and documentation of endoscopic and microscopic procedures.

5.4 Contraindications

No contraindications relating directly to the medical device are currently known. The responsible physician must decide whether the foreseen application is admissible based on the general condition of the patient.

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5.5 Patient profile

Use of the IMAGE1 S™ is not limited to a certain patient profile (sex, age, weight etc.). The attending physician must assess a patient's health with regard to the endoscopic intervention and declare it adequate.

5.6 User profile (specialist physician)

- Recognized medical skills for endoscopic interventions in the relevant specialization
- Adequate powers of comprehension to assess the situation rationally at all times during the application
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- No physical impairments that could diminish perception of activation and visual message

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
5.7 User profile (assistant)

- Completed training as medical assistant
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- No physical impairments that could diminish perception of activation and visual message

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5.8 Safety precautions at the site of installation

The unit may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations. The system is not intended for use in areas where there is a risk of explosion. This means, for example:

Wherever easily combustible and explosive inhalation anesthetics and their mixtures are used, the device may not be operated in the demarcated hazard zone. This also applies to easily combustible and potentially explosive chemicals, e.g., skin disinfectants and fast-acting surface disinfectants.

The unit is equipped with a connector for potential equalization. It should be connected in accordance with the applicable national regulations.

5.9 Safety precautions when operating the device

It is the user's responsibility to make sure the device is safe and operating properly before using the device.

During treatment with the equipment, the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

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6.2 Unpacking the equipment

Carefully remove the device and accessories from the packaging. Check for missing items and evidence of damage during shipment. Please file any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these may be useful if the unit has to be transported.

6.3 Preparing the device for use

6.3.1 Installation

- NOTE:** Up to three IMAGE1 S[™] LINK modules can be connected to the IMAGE1 S CONNECT[®].
- NOTE:** Do not cover the lateral air vents. Maintain a minimum distance of 6 cm to other devices.

Place the IMAGE1 S CONNECT[®] and the corresponding IMAGE1 S[™] LINK modules on a flat surface.

6.3.2 Connecting the potential equalization

- NOTE:** The unit's ground line should be installed by a qualified electrician.

The IMAGE1 S CONNECT[®] and the IMAGE1 S[™] LINK modules are equipped with connectors for potential equalization (see fig. A and B).

1. Connect the potential equalization line to the connectors for potential equalization on all modules.
2. Connect the potential equalization line to the connector in the treatment room.

6.3.3 Connecting the IMAGE1 S[™] LINK module with IMAGE1 S CONNECT[®]

Using the LINK cable supplied, connect the LINK module to the IMAGE1 S CONNECT[®] (see fig. C).

- NOTE:** The IMAGE1 S CONNECT[®] inputs are numbered. Please connect the LINK modules in ascending order.

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
6.3.4 Connecting SCB

Pull back the protection device of the SCB connector and insert the connector into one of the SCB sockets of the TC200 IMAGE1 S CONNECT® or into the SCB socket of the TC201 IMAGE1 S CONNECT® II. Connect the other end of the cable to the KARL STORZ-SCB (KARL STORZ Communication Bus) control unit or other SCB units (see KARL STORZ SCB control NEO System Instruction Manual) (see fig. D).

6.3.5 Establishing a KS HIVE connection

Plug the OR1™ patch cable into the LAN interface of the TC200 IMAGE1 S CONNECT® or into the LAN interface of the TC201 IMAGE1 S CONNECT® II (see fig. E or F). Connect the other end of the cable with another KS HIVE-compatible device.*

* At present only one of the following devices can be used:

- TL400 Power LED Rubina
- UR500B, UR500C ARTIP CRUISE™

6.3.6 Connecting the monitor

Connect the monitor with the IMAGE1 S CONNECT® or CONNECT II (see fig. G). Depending on the types of connectors provided on the monitor, this connection can be generated using either:

- DVI-D or 3G-SDI output for IMAGE1 S CONNECT® TC200
- DisplayPort, 12G-/3G-SDI or DVI-D output for IMAGE1 S CONNECT® II TC201

The primary working screen should be connected to DVI1 for TC200.

The primary working screen should be connected to DP1 for TC201.

- ❗ **NOTE:** The connected monitor must support a resolution of 1920 x 1080 pixels.
- ❗ **NOTE:** The connected monitor must support a resolution of 3840 x 2160 pixels to support 4K UHD.
- ❗ **NOTE:** A compatible 3D monitor must be connected when using a D3-LINK in conjunction with 3D video endoscopes (see chapter 11.11).

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

6.3.7 Connecting USB devices

The USB ports of the IMAGE1 S CONNECT® are suitable for connecting up the following peripherals (see fig. H):

- Keyboard/Mouse
- Printer
- Footswitch
- USB to ACC adaptor
- USB flash drives, FAT32
- IMAGE1 PILOT TC014

6.3.8 Connecting the power supply

Before connecting the devices to the electrical supply line, check that the values on the power sockets concur with the values on the devices' identification plates.

1. Insert the device plug of the power cord into the power supply cord receptacle of the IMAGE1 S CONNECT® (see fig. I).
2. Insert the device plug of the power cord into the power supply cord receptacle of the IMAGE1 S™ LINK modules (see fig. J).
3. Insert the modules' power cords into the power sockets.

WARNING: Always use a hospital grade power cord with this camera.

6.3.9 Switching on the IMAGE1 S CONNECT®

Switch on IMAGE1 S CONNECT® (see fig. K).

NOTE: The connected LINK modules are turned on automatically.

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6.3.10 Connecting USB storage media

Connect the USB storage media with the IMAGE1 S CONNECT® to store still images and/or videos (see fig. L).

NOTE: The USB storage medium must have a FAT 32 format.

6.3.11 Connecting the camera head/ video endoscope

Insert the camera head connector into the receptacle of the relevant LINK module (see fig. M).

NOTE: Please verify the values set for the primary and secondary sources as well as the layout beforehand (see p. 23).

NOTE: Please note the compatibility list of the respective LINK modules (see chap. 11).

WARNING: Always ensure the connector is completely clean and free of debris. If moisture is present, dry thoroughly with sterile towel prior to insertion.

6.3.12 Performing the white balance

NOTE: Perform a white balance before every use, in order to adjust the camera's perception of white and the ability to display colors properly.

1. Ensure that an imaging instrument* is connected to an IMAGE1 S™ LINK module and showing a live image on the monitor.
2. Connect the light source to the imaging instrument via the light cable.
3. Point the tip of the imaging instrument at any clean white surface (e.g. white gauze).
4. Ensure that the light source is emitting enough light to sufficiently illuminate the white surface. The monitor should display an entire white and well illuminated live image.

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5. Perform the white balance by pressing the  button on the front of the IMAGE1 S CONNECT® or by pressing the  in the Live Menu.
6. Look at the monitor to determine whether the white balance was successful or whether it failed.

*A camera head with an attached endoscope, video endoscope, or exoscope

6.3.13 Test for proper functioning

1. Inspect the IMAGE1 S™ module and the camera head and/or the video endoscope for external damage.
2. Ensure that the camera head cable and/or the video endoscope cable is not damaged or kinked.
3. Point the camera head and/or the video endoscope towards an object and check the quality of the display on the monitor.
4. When using 3D video endoscopes, check to make sure the monitor's display is correct. Please ensure that you use the correct polarized eyewear for the monitor.

 **NOTE:** The device and accessories may only be used once the test for proper functioning has been successfully completed.

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- USB-Stick
- IMAGE1 PILOT (3D-Maus)
- USB-Hub
- 3D-Fußschalter

* KS HIVE connection is also possible for ARTIP CRUISE™

** USB connection is also possible for

- FLASH memory stick
- IMAGE1 PILOT (3D-Mouse)
- USB-Hub
- 3D-Footswitch

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

7 Operating instructions

7.1 Basic information on operation

The use of IMAGE1 S™ is divided into three menu stages:

- Dashboard (see 7.2)
- Live menu (see 7.3)
- Setup menu (see 7.4)

7.1.1 Keyboard and mouse

The keyboard and mouse can be used to navigate the menus. A keyboard is required to enter patient data.

7.1.2 Camera head buttons

All IMAGE1 camera heads feature three buttons for retrieving programmed functions and for menu control (see fig. N):

- ⑲ scroll up: Scrolls up the menus, options.
- ⑳ scroll down: Scrolls down the menus, options.
- ㉑ menu/select: Displays the live menu. Once this menu is displayed, this button activates the highlighted option (identical to ENTER on a computer keyboard).

7.2 Dashboard

The Dashboard displays the devices, data storage devices and LINK modules connected to IMAGE1 S CONNECT® (see fig. O). Unconnected devices are shown in gray.

- NOTE:** The Dashboard is shown when no camera head/video endoscope is connected.
- NOTE:** Depending on the configuration of your system, individual icons may be deactivated or unavailable. The configuration can be modified in the Setup Menu.
- NOTE:** The icons display the current settings.

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7.3 Live Menu

The live menu allows you to choose between different camera functions and change the settings of connected devices (see fig. P). Unconnected devices are shown in gray.

You can access the live menu by pressing the enter button on the keyboard or by pressing the menu button  on the camera head/ video endoscope or IMAGE1 PILOT (TC014) left side button.

-  **NOTE:** The icons in the live menu are displayed as soon as a camera head is connected.
-  **NOTE:** Depending on the configuration of your system, individual icons may be deactivated or unavailable. The configuration can be modified in the Setup Menu.
-  **NOTE:** The icons display the current settings.

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7.3.1 Description of functions of the symbols used

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
7.5 Patient Data Privacy

NOTE: By default the IMAGE1 S™ system does not contain any patient data. The patient data management functions (e.g. enter, access, display, store, or transmit) have to be enabled by the user.

NOTE: By enabling the patient data management functions the user accepts the responsibility for all patient related data managed with the help of IMAGE1 S™ according to applicable national data protection laws.

NOTE: Do not record sensitive patient information including patient's face, name, social security number, etc. with IMAGE1 S™ camera system.

NOTE: KARL STORZ assumes no legal liability or responsibility for loss, correctness, damage or exposure of patient data managed with the help of IMAGE1 S™.

In the Patient Data Privacy settings, you can configure all settings related to patient information. Patient information of any kind is only available if the corresponding settings have been configured.

Configuration of the **Patient Data Privacy** menu is password-protected as standard.

When activating it for the first time, an individual password must be assigned (see fig. Q and R). The password must meet the following conditions:

- at least eight characters
- of which at least one capital letter and one lower case letter
- at least one digit and at least one special character.

After entering the password once, it must be confirmed in the second line. If both lines are valid and the password is confirmed, the menu can be edited.

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The Increased Patient Data Privacy item provides you with preset cyber security settings. When this setting is enabled (= active), all Patient Data Privacy sub configurations are disabled and set to the default settings. The affected sub configurations and their default values are:

- ▷ Display Duration
 - You can determine how long the patient data is displayed in the Patient Region of the screen.
 - Default: Never
- ▷ Show Patient Tiles
 - You can determine whether the Patient Tile will be shown on the Dashboard and Live Menu or not.
 - Default: No
- ▷ Enable the printing of Patient Data on Prints
 - You can determine whether prints will be labeled with patient information or not.
 - Default: No
- ▷ Add Patient Names to Files Captured
 - You can determine whether media files stored will be labeled with the patient information.
 - Default: No
- ▷ Patient Data Storage
 - You can specify how long data is stored on the device.
 - Default: 1 day
- ▷ View/Save Logged Events
 - You can view the Audit log by activating the Open Log button. Audit Logs can be exported via the 'Save As...' function.

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The Submit to Service / Reset Password option provides access to the settings of the Patient Data Privacy menu if the current password is no longer available.

Please note: For security reasons all patient data will be deleted from the device and the device will be reset to the default patient data privacy settings.

A new password must be set before patient data can be entered. Use of this function is also recommended if the device is to be sent in for maintenance work.

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8 Reprocessing

8.1 General safety information

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

CAUTION: When preparing and using solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: Do not allow water or any liquid to drip inside the unit(s).

CAUTION: Damage to product. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be requested from hygiene@karlstorz.com. *Not available for U.S. Customers

CAUTION: National laws and regulations must be observed.

NOTE: The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

NOTE FOR US CUSTOMERS: The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' is not released for use in the United States.

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8.2 Wipe-down disinfection of devices

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry low lint cloth.

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9.1 Maintenance and safety check

9.1.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

9.1.2 Safety check

Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, the following checks must be carried out on this device at least once a year by persons who, due to their training, knowledge and practical experience, are able to properly carry out such safety checks and who are not subject to any instructions as far as such checking activities are concerned.

Visual inspection:

- Check the device and accessories for any mechanical damage which may impair functionality.
- Check that the labels relevant to safety are legible.

Electric measurements:

- Protective grounding measured according to IEC 62353 (current edition), and meets the limits specified in the standard.
- Touch Current, Earth Leakage Current, and Patient Leakage currents, measured according to IEC 62353 (current edition), and meets the limits specified in the standard.
- Testing at system level in the cart includes patient leakage testing. The IEC 62353 patient leakage measurement includes testing of special (fault) conditions of mains on Applied Parts, and is conducted with one Type CF applied part measured in turn as per Clause 5.3.4.3.1

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& Figure B2 of the IEC 62353, (limit 50 µA). When tested according to the IEC 60601-1 3.1 ed. at system level in the cart, it includes total leakage current of testing in normal conditions (limit 50 µA) as per clause # 16.6.3

Test for proper functioning:

- A functional check must be carried out in accordance with the instruction manual (see § 6.3.13).
- The safety check must be entered in the unit log-book and the results of the check must be documented.
- If the device does not function reliably and/or safely, it must be repaired.

9.2 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by KARL STORZ: all repair work must employ genuine KARL STORZ parts only.

9.3 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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9.4 Repairing devices

Only personnel authorized by the manufacturer may repair devices. Interventions described in these instructions for use are exempt from this rule.

In Germany, you can contact the following address for repairs:

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Take-Off Gewerbepark 83
78579 Neuhausen, Germany

Service hotline: +49 7461 708-980
E-mail: technicalsupport@karlstorz.com

In other countries, the respective KARL STORZ branches or specialist dealers are responsible.

9.5 Infection prevention

Contaminated devices may not be shipped. To prevent contact infections and airborne infections, devices must first be decontaminated. KARL STORZ reserves the right to send back contaminated devices.

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9.6 Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- the electrical installation of the location in which the device is connected and operated conforms to the applicable statutes and standards, and
- the unit has been used in accordance with its operating instructions at all times.

9.7 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see 'Subsidiaries' section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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10.2 Technical documentation

On request, the manufacturer will provide available parts lists, descriptions, instructions for adjustments, and other documentation deemed useful for suitably qualified user personnel authorized by the manufacturer when repairing device components that have been designated as repairable by the manufacturer.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

10.3 Standard compliance

According to IEC 60601-1 + A1 + A2, IEC 60601-2-18, UL 60601-1, CSA C22.2 No. 601.1-M90:

- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part of type CF defibrillation-proof 
- Type of protection against moisture: IPX 0

According to IEC 60601-1-2:

Please observe the information on electromagnetic compatibility in § 12 (p. 66).

10.4 Directive compliance for IMAGE1 S™ (TC200, TC201, TC300, TC301, TC302 and TC304)

According to Medical Device Directive (MDD): This medical device belongs to Class I

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

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10.5 Software ownership and licensing

With respect to KARL STORZ products containing software components, Customer has a non-exclusive, limited, non-transferable license to use the programmed logic, computer programs and/or software supplied by KARL STORZ, in connection with, and incorporated into, KARL STORZ Products System (Software) internally, only in the form in which delivered to Customer and for the sole purpose of operating in accordance with KARL STORZ's written instructions for such products sold to Customer (and for no other product). The Software, and all modifications, enhancements and upgrades thereto, will, at all times, remain the property of KSEA. Customer may not duplicate, copy, reverse-engineer, de-compile, or disassemble the Software or in any way modify the Software. Customer has no right to, and may not, create derivatives of the software, and Customer may not attempt to copy, create or re-create the source code of the Software. Any and all such modifications or enhancements to the Software by Customer, in contravention of this license, will immediately become the sole property of KSEA.

Customer hereby acknowledges and agrees (I) that the purchase of KARL STORZ products does not constitute a sale of the Software, (II) that the Software is the property of KARL STORZ, (III) that Customer neither owns nor acquires any claim or intellectual property right in or to the Software, and (IV) KARL STORZ retains all right, title, and interest in and to the Software, at all times, regardless of the form or media in or on which the original or other copies of the Software may exist. In the event of a failure of Customer, or its agents, employees or representatives, to comply with any terms and conditions of the License herein granted, the License will, without any further action by KARL STORZ or any other party, immediately terminate.

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Attribution Notice

This device contains open source software. You can view a list of the open source software included in the product, and the open source software license terms, on the IMAGE1 S CONNECT® device itself. From the Dashboard, select 'Setup/General', then 'Software Licenses'.

To obtain a copy of the source code being made publicly available by KARL STORZ Imaging, Inc. ('KARL STORZ') related to open source software used in this KARL STORZ product, you should send your request, along with product part number and serial number, in writing by email to opensource@karlstorz.com or by postal mail to:

KARL STORZ Imaging, Inc.
Attention: Open Source Software Inquiries
1 S. Los Cameros
Goleta, CA 93117
USA

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**11 Spare parts,
recommended accessories**

**11.1 Compatible camera heads for
IMAGE1 S™ H3-LINK TC300
with S technologies**

Name	Article number
IMAGE1 S™ H3-Z Three- Chip FULL HD Camera Head	TH100
IMAGE1 S™ H3-ZI Three- Chip FULL HD Camera Head	TH101
IMAGE1 S™ H3-Z FI Three- Chip FULL HD Camera Head	TH102
IMAGE1 S™ H3-P Three- Chip FULL HD Pendulum Camera Head	TH103
IMAGE1 S™ H3-ZA Three- Chip FULL HD Camera Head	TH104
IMAGE1 S™ H3-M COVIEW® Three-Chip FULL HD Camera Head	TH106

**11.2 Compatible camera heads for
IMAGE1 S™ H3-LINK TC300
without S technologies**

Name	Article number
IMAGE1 H3-Z Three-Chip FULL HD Camera Head	22220055-3
IMAGE1 H3-ZI Three-Chip FULL HD Inline Camera Head	22220053-3
IMAGE1 H3-Z FI Three-Chip ICG FULL HD Camera Head	22220085-3*
IMAGE1 H3-P Three-Chip FULL HD Pendulum Camera Head	22220056-3
IMAGE1 H3-ZA Three-Chip FULL HD Camera Head	22220061-3
IMAGE1 H3-FA Three-Chip FULL HD Camera Head	22220060-3

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Information on electromagnetic compatibility (EMC)

12 Information on electromagnetic compatibility (EMC)

The IMAGE1 S™ is suitable for use within close proximity to an active HF electro-surgical device in professional healthcare facility environment. Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

WARNING: This equipment is MR Unsafe. Keep the device away from the Magnetic Resonance Imaging (MRI) Scanner Room and mobile MRI scanner devices.

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

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System Cables and Maximum Lengths used for EMC Compliance				
Type	Shielded	Max Length [m]	Ferrite	Used for
LINK cable	Yes	2	No	Interconnection of inter-module LINK signals
BNC to BNC	Yes	10	No	Connection of 3G/12G-SDI video signals
DVI	Yes	3	No	Connection of DVI signals
DisplayPort (DP)	Yes	3	No	Connection of DisplayPort signals
6-pin mini DIN	Yes	1	No	Connection of SCB signals
USB	Yes	2	No	Connection to printer, storage device, keyboard
Power cord	No	3	No	Connection of Device to AC Mains
Ethernet	Yes	2	No	Connection to LAN (RJ-45) port
IMAGE1 PILOT TC014	Yes	5	No	Connection to VITOM® 3D
External synchronization cable	Yes	1	No	Connection to compatible light source

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Information on electromagnetic compatibility (EMC)

WARNING: The use of an accessory or cable with the IMAGE1 S™ other than those specified in this manual may result in increased emissions or decreased immunity of the IMAGE1 S™. When using an accessory or cable with IMAGE1 S™ other than those specified in this manual, it becomes the responsibility of the user of the IMAGE1 S™ to determine compliance with IEC 60601-1-2 when using this item.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the IMAGE1 S™, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

NOTE: The emission characteristics of this equipment make it suitable for use in industrial areas as well as in hospitals (CISPR 11 Class A) and other professional healthcare environment. If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio communication service. The user might need to take mitigation measures, such as relocating or re-orienting equipment.

NOTE: HF surgical instrument compatibility

This equipment has been evaluated for compatibility with high-frequency (HF) surgical equipment. It has been tested according to Annex BB of the standard IEC 60601-2-2.

NOTE: RFID compatibility

This equipment has been evaluated for compatibility with Radio-Frequency Identification (RFID). It has been tested according to the standard AIM 7351731.

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Table 1 – Compliance level for Immunity tests
Guidance and manufacturer's declaration – electromagnetic immunity

The IMAGE1 S™ is intended for use in the electromagnetic environment specified below. The user of IMAGE1 S™ should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV Contact discharge ± 15 kV Air discharge	± 8 kV Contact discharge ± 15 kV Air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1 kV ± 1 kV for input and output lines 100 kHz repetition	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage Dip:</u> Dip to 0 % for 1 cycle @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	<u>Voltage Dip:</u> Dip to 0 % for 1 cycle @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the IMAGE1 S™ require continued operation during power mains interruptions, it is recommended that the IMAGE1 S™ be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50 Hz/60 Hz	30 A/m at 50 Hz/60 Hz	If image distortion occurs, it may be necessary to position the IMAGE1 S™ further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3 V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3 V _{rms} on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V _{rms} in ISM bands	3 V _{rms} on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V _{rms} in ISM bands	

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Table 2

Test levels for Proximity fields from RF wireless communications equipment

Test Frequency MHz	Band MHz	Service	Modulation	Immunity Test Level V/m	Compliance level V/m
385	380 – 390	TETRA 400	Pulse modulation 18 Hz	27	27
450	430 – 470	GMRS 460 FRS 460	FM $\pm$ 5 kHz deviation 1 kHz Sine Wave	28	28
710	704 – 787	LTE Band 13 & 17	Pulse modulation 217 Hz	9	9
745					
780					
810	800 – 960	GSM 800/900 TETRA 800 iDEN 820 CDMA 850 LTE Band 5	Pulse modulation 18 Hz	28	28
870					
930					
1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	28	28
1845					
1970					
2450	2400 – 2570	Bluetooth WLAN 802.11 b/g/n RFID 2450 LTE Band 7	Pulse modulation 217 Hz	28	28
5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9	9
5500					
5785					

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Table 3 – Test levels for Radiated and conducted Immunity Tests

Guidance and manufacturer's declaration – electromagnetic immunity

The IMAGE1 S™ is intended for use in the electromagnetic environment specified below. The user of IMAGE1 S™ should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} on 150 kHz to 80 MHz 6 V _{rms} in ISM bands between 0.15 MHz and 80 MHz	3 V _{rms} on 150 kHz to 80 MHz 6 V _{rms} in ISM bands between 0.15 MHz and 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the IMAGE1 S™, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b .
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m 80 MHz to 2.7 GHz	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.7 GHz Interference may occur in the vicinity of equipment marked with the following symbol: 

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects, and persons.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the IMAGE1 S™ is used exceeds the applicable RF compliance level above, the IMAGE1 S™ should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the IMAGE1 S™.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 – Emission class and group

Guidelines and manufacturer's declaration – electromagnetic emissions

The IMAGE1 S™ is intended for use in the electromagnetic environment specified below. The customer or user of the IMAGE1 S™ should ensure that it is used in such an environment.

RF emissions as per CISPR 11	Group 1	The IMAGE1 S™ uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class A	The IMAGE1 S™ is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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Table 5

Recommended separation distances between portable and mobile RF communications equipment and the IMAGE1 S™

The IMAGE1 S™ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the IMAGE1 S™ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the IMAGE1 S™ as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.

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13 Appendix

13.1 Compatible Veterinary video endoscopes for IMAGE1 S™ X-LINK TC301 with S technologies

Name	Article number (PAL/NTSC)
Specialty VET	
Veterinary Video Endoscope	60914NKS/PKS
Veterinary Video Endoscope	60714NKS/PKS
Veterinary Video Endoscope	60511NKS/PKS
Veterinary Video Endoscope	60332NKS/PKS
Veterinary Video Endoscope	60130NKS/PKS
Veterinary Video Endoscope	60125NKS/PKS
Veterinary Video Endoscope	60118NKS/PKS
Veterinary Video Endoscope	60278VS/VSU

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Carneros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N.W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorz@cuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Av. Ejercito Nacional No. 453 Piso 2,
Colonia Granada, Alcaldia Miguel Hidalgo
C.P. 11520 Ciudad de México
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA – Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Stamveien 1
1483 Hagan, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
127 39 Skärholmen, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se
KARL STORZ Endoscopy Suomi Oy
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestučio st. 59 / Lenkioji st. 27
08124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 67 000
E-Mail: info-lt-lv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk
KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopia Nederland B. V.
Displayweg 2
3821 BT Amerfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com
KARL STORZ Endoscopia France S.A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.com

KARL STORZ Endoscopia Ibérica S.A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-it@karlstorz.com

KARL STORZ Croatia d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info-hrv@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Cesta v Gorice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Polska Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy China Ltd.
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch
Room 701A Building 5 & Room 501 Building 7,
No. 3000 Longdong Avenue, Pilot Free Trade Zone,
201203, Shanghai, People's Republic of China
Phone: +86 21 60339888, Fax: +86 21 60339808
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
Room 803-805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch
Room 02B & 03 & 04A, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 87321281, Fax: +86 20 87321286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Indonesia
Sinarmas MSG Tower Level 37
Jl. Jend. Sudirman No. Kav. 21
Jakarta Selatan
DKI Jakarta 12920
E-Mail: infoindonesia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info-jp@karlstorz.com

KARL STORZ Endoscopy New Zealand Ltd.
31 Morningside Drive Mt Albert
Auckland, 1025, New Zealand
PO Box 56 511, Dominion Rd
Auckland, 1446, New Zealand
Phone: +64 9 846 6044, Fax: +64 9 846 6808
Toll free: +64 508 84 84 84 (New Zealand only)
E-Mail: sales@karlstorz.co.nz

KARL STORZ Endoscopy Australia Pty. Ltd.
68 Waterloo Road, Macquarie Park NSW 2113
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: karlstorz@karlstorz.com.au
www.karlstorz.com

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