Abbott Rapid Diagnostics Jena GmbH
Orlaweg 1
07743 Jena
Germany

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EN

Introduction

The Coronavirus disease (COVID-19) is an infectious disease caused by a newly discovered coronavirus, severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2)¹. The SARS-CoV-2 is a β -coronavirus, which is an enveloped non-segmented positive-sense RNA virus². It is spread by human-to-human transmission via droplets or direct contact, and infection has been estimated to have a mean incubation period of 6.4 days and a basic reproduction number of 2.24-3.58. Among patients with pneumonia caused by SARS-CoV-2, fever was the most common symptom, followed by cough³. The main IVD assays used for COVID-19 employ real-time reverse transcriptase-polymerase chain reaction (RT-PCR) that takes a few hours⁴. The availability of a cost-effective, rapid point-of-care diagnostic test is critical to enable healthcare professionals to aid in the diagnosis of patients and prevent further spread of the virus⁵. Antigen tests will play a critical role in the fight against COVID-19⁶.

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Intended Use

Panbio™ COVID-19 Ag Rapid Test Device is an *in vitro* diagnostic rapid test for the qualitative detection of SARS-CoV-2 antigen (Ag) in human nasopharyngeal swab specimens from individuals who meet COVID-19 clinical and / or epidemiological criteria. Panbio™ COVID-19 Ag Rapid Test Device is for professional use only and is intended to be used as an aid in the diagnosis of SARS-CoV-2 infection. The product may be used in any laboratory and non-laboratory environment that meets the requirements specified in the Instructions for Use and local regulation. The test provides preliminary test results. Negative results don't preclude SARS-CoV-2 infection and they cannot be used as the sole basis for treatment or other management decisions. Negative results must be

ENGLISH

combined with clinical observations, patient history, and epidemiological information. The test is not intended to be used as a donor screening test for SARS-CoV-2.

Kit Variants

- 41FK10 No 2D barcode printed on the contained test devices
- 41FK20 Contains test devices with a 2D barcode printed on the test device, which encodes traceability information for the product

Materials Provided

- 25 Test devices with desiccant in individual foil pouch
- Buffer (1 x 9 ml/bottle)
- 25 Extraction tubes
- 25 Extraction tube caps
- 1 Positive control swab
- 1 Negative control swab
- 25 Sterilized nasopharyngeal swabs for sample collection
- 1 Tube rack
- 1 Quick Reference Guide (Nasopharyngeal)
- 1 Instructions for use

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Storage and Stability

1. The test kit should be stored at a temperature between 2-30 °C. Do not freeze the kit or its components.

Note: When stored in a refrigerator, all kit components must be brought to room temperature (15-30 °C) for a minimum of 30 minutes prior to performing the test. Do not open the pouch while components come to room temperature.

2. The Buffer bottle may be opened and resealed for each assay. The Buffer cap should be firmly sealed between each use. The Buffer is stable until expiration date if kept at 2-30 °C.
3. Perform the test immediately after removing the test device from the foil pouch.
4. Do not use the test kit beyond its expiration date.
5. The shelf life of the kit is as indicated on the outer package.
6. Do not use the test kit if the pouch is damaged or the seal is broken.
7. Direct swab specimens should be tested immediately after collection. If immediate testing is not possible, the swab specimen can be kept in an extraction tube filled with extraction buffer (300 µl) at room temperature (15-30 °C) for up to two hours prior to testing.

Warnings

1. For *in vitro* diagnostic use only. Do not reuse the test device and kit components.
2. These instructions must be strictly followed by a trained healthcare professional to achieve accurate results. All users have to read the instruction prior to performing a test.
3. Do not eat or smoke while handling specimens.
4. Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
5. Avoid splashing or aerosol formation of specimen and buffer.
6. Clean up spills thoroughly using an appropriate disinfectant.
7. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials (i.e. swab, extraction tube, test device) in a biohazard container as if they were infectious waste and dispose according to applicable local regulations.
8. Do not mix or interchange different specimens.
9. Do not mix reagent of different lots or those for other products.
10. Do not store the test kit in direct sunlight.
11. To avoid contamination, do not touch the head of provided swab when opening the swab pouch.
12. The provided sterilized swabs in the package should be used only for nasopharyngeal specimen collection.
13. To avoid cross-contamination, do not reuse the sterilized swabs for specimen collection.
14. Do not dilute the collected swab with any solution except for the provided extraction buffer.
15. The buffer contains <0.1% sodium azide as a preservative which may be toxic if ingested. When disposed of through a sink, flush with a large volume of water.⁷

ENGLISH

16. Do not use the positive or negative control swab for specimen collection.

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Specimen Collection & Extraction

1. Tilt the patient's head back slightly about 45°-70° to straighten the passage from the front of the nose.
2. Insert the swab with a flexible shaft through the nostril parallel to the palate.
⚠ Caution: Use dedicated nasopharyngeal swab for specimen collection.
3. Swab should reach depth equal to distance from nostrils to outer opening of the ear.
⚠ Caution: If resistance is encountered during insertion of the swab, remove it and attempt insertion in the opposite nostril.
4. Gently rub and roll the swab, 3-4 times. Leave the swab in place for several seconds to absorb secretions.
5. Slowly remove swab while rotating it and insert into the extraction tube.
6. Swirl the swab tip in the buffer fluid inside the extraction tube, pushing into the wall of the extraction tube at least five times and then squeeze out the swab by squeezing the extraction tube with your fingers.
7. Break the swab at the breakpoint and close the cap of extraction tube.

Reaction with Test Device

1. Open the dropping nozzle cap at the bottom of the extraction tube.
2. Dispense 5 drops of extracted specimens vertically into the specimen

ENGLISH

well (S) on the device. Do not handle or move the test device until the test is complete and ready for reading.

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1. **Negative result:** The presence of only the control line (C) and no test line (T) within the result window indicates a negative result.
2. **Positive result:** The presence of the test line (T) and the control line (C) within the result window, regardless of which line appears first, indicates a positive result.

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3. **Invalid result:** If the control line (C) is not visible within the result window after performing the test, the result is considered invalid.

Test Limitations

1. The contents of this kit are to be used for the professional and qualitative detection of SARS-CoV-2 antigen from nasopharyngeal swab. Other specimen types may lead to incorrect results and must not be used.
2. Failure to follow the instructions for test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
3. A negative test result may occur if the specimen was collected, extracted or transported improperly. A negative test result does not eliminate the possibility of SARS-CoV-2 infection and should be confirmed by viral culture or a molecular assay.
4. Positive test results do not rule out co-infections with other pathogens.
5. Test results must be evaluated in conjunction with other clinical data available to the physician.
6. Reading the test results earlier than 15 minutes or later than 20 minutes may give incorrect results.
7. Panbio™ COVID-19 Ag Rapid Test Device is not intended to detect from defective (non-infectious) virus during the later stages of viral shedding that might be detected by PCR molecular tests.⁸
8. Positive results may occur in cases of infection with SARS-CoV.

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2. External Quality Control:

The controls are specifically formulated and manufactured to ensure performance of the Panbio™ COVID-19 Ag Rapid Test Device and are used to verify the user's ability to properly perform the test and interpret the results. The Positive Control contains recombinant

ENGLISH

SARS-CoV-2 nucleocapsid protein, which is not contagious. The Positive Control will produce a positive test result and has been manufactured to produce a visible test line (T). The Negative Control will produce a negative test result. Control swabs are not specific for a particular Panbio™ COVID-19 Ag Rapid Test Device lot and may be used between test device lots until the swabs' expiry dates. Good laboratory practice suggests the use of positive and negative controls to ensure that:

- Test reagents are working, and
 - The test is correctly performed.
- The external controls can be run under any of the following circumstances:
- By a new operator prior to performing testing on patient specimens,
 - When receiving a new test shipment,
 - At periodic intervals as dictated by local requirements, and/or by the user's Quality Control procedures.

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Clinical performance of Panbio™ COVID-19 Ag Rapid Test Device was determined by testing 140 positive and 445 negative specimens for SARS-CoV-2 antigen (Ag) to have a sensitivity of 91.4% (95% CI: 85.5-95.5%) and specificity of 99.8% (95% CI: 98.8-100%). Clinical specimens were determined to be positive or negative using an FDA EUA RT-PCR reference method.

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		PCR Test Result		
		Positive	Negative	Total
Panbio™ COVID-19 Ag Rapid Test Device Results	Positive	128	1	129
	Negative	12	444	456
	Total	140	445	585
		Sensitivity	Specificity	Overall Percent Agreement
		91.4%	99.8%	97.8%
		[85.5%; 95.5%]	[98.8%; 100%]	[96.2%; 98.8%]

- Performance data was calculated from a study of individuals suspected of exposure to COVID-19 or who have presented with symptoms in the last 7 days.

- Stratification of the positive specimens post onset of symptoms or suspected exposure between 0-3 days has a sensitivity of 94.9% (95% CI: 82.7-99.4%; n=39) and 4-7 days has a sensitivity of 90.1% (95% CI 82.5-95.1%; n=101).
- Positive agreement of the Panbio™ COVID-19 Ag Rapid Test Device is higher with samples of Ct values ≤30 with a sensitivity of 97.6% (95% CI: 93.2-99.5%) and Ct values ≤33 with a sensitivity of 94.1% (95% CI: 88.7-97.4%). As suggested in References 8-10, patients with Ct value >30 are no longer contagious.^{8, 9, 10}

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Clinical performance of Panbio™ COVID-19 Ag Rapid Test Device was determined by testing 483 asymptomatic subjects for SARS-CoV-2 antigen (Ag). Clinical specimens were determined to be positive or negative using an FDA EUA RT-PCR reference method. The positive results (n=50) were stratified by the comparator method cycle threshold (Ct) counts and assessed to better understand the correlation of product performance, as a surrogate for the amount of virus present in the clinical sample. A lower Ct value corresponds to a higher virus concentration. As presented in the table below, the positive agreement increases with lower Ct values.

This clinical data was collected using a nasal swab specimen and is presented for information purposes only to reflect product performance in asymptomatic subjects. This product kit contains only accessories for nasopharyngeal sample collection.

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As suggested in References 8-10, patients with Ct value >30 are no longer contagious.^{8, 9, 10}

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4. Hook Effect

There is no hook effect at $1.0 \times 10^{5.8}$ TCID₅₀/ml of SARS-CoV-2

ENGLISH

which was isolated from a COVID-19 confirmed patient in Korea.

5. Cross Reactivity

Cross-reactivity of Panbio™ COVID-19 Ag Rapid Test Device was evaluated by testing 27 viruses and 21 other microorganisms. The final test concentrations of viruses and other microorganisms are documented in the Table below. The following viruses and other microorganisms except the Human SARS-coronavirus Nucleoprotein have no effect on the test results of Panbio™ COVID-19 Ag Rapid Test Device. Panbio™ COVID-19 Ag Rapid Test Device has cross-reactivity with Human SARS-coronavirus Nucleoprotein at a concentration of 25 ng/ml or more because SARS-CoV has high homology (79.6%) to the SARS-CoV-2.

No.	Types of Specimen	Cross Reaction Substance	Final Test Concentration	Test Result
1	Virus	Adenovirus Type3	2.0 X 10 ^{6.5} TCID ₅₀ /ml	No cross reaction
2		Adenovirus Type7	2.0 X 10 ^{4.75} TCID ₅₀ /ml	No cross reaction
3		Echovirus2	1.0 X 10 ^{6.5} TCID ₅₀ /ml	No cross reaction
4		Echovirus11	2.0 X 10 ^{5.25} TCID ₅₀ /ml	No cross reaction
5		Human herpesvirus (HSV) 1	2.0 X 10 ^{6.25} TCID ₅₀ /ml	No cross reaction
6		Human herpesvirus (HSV) 2	2.0 X 10 ^{4.75} TCID ₅₀ /ml	No cross reaction
7		Mumps Virus Ag	2.0 X 10 ^{3.5} TCID ₅₀ /ml	No cross reaction
8		Influenza virus A (H1N1) Strain (A/Virginia/ATCC1/2009)	2.6 X 10 ^{5.0} PFU/ml	No cross reaction
9		Influenza virus A (H1N1) Strain (A/WS/33)	5.0 X 10 ^{7.25} TCID ₅₀ /ml	No cross reaction
10		Influenza virus A(H3N2) Strain (A/Hong Kong/8/68)	N/A*	No cross reaction
11		Influenza virus B Strain (B/Lee/40)	2.0 X 10 ^{5.25} TCID ₅₀ /ml	No cross reaction
12-14		Parainfluenza Type 1, Parainfluenza Type 2, Parainfluenza Type 3	N/A*	No cross reaction

ENGLISH

No.	Types of Specimen	Cross Reaction Substance	Final Test Concentration	Test Result
15	Virus	Parainfluenza Type 4A	$1.97 \times 10^{7.0}$ PFU/ml	No cross reaction
16		Respiratory syncytial virus (RSV) type A	$4.22 \times 10^{5.0}$ TCID ₅₀ /ml	No cross reaction
17		Respiratory syncytial virus (RSV) type B	$5.62 \times 10^{5.0}$ TCID ₅₀ /ml	No cross reaction
18		HCoV-HKU1	10 µg/ml	No cross reaction
19		Rhinovirus A16	$8.8 \times 10^{5.0}$ PFU/ml	No cross reaction
20		HCoV-NL63	$1.7 \times 10^{5.0}$ TCID ₅₀ /ml	No cross reaction
21		HCoV-OC43	$8.9 \times 10^{5.0}$ TCID ₅₀ /ml	No cross reaction
22		HCoV-229E	$1.51 \times 10^{6.0}$ TCID ₅₀ /ml	No cross reaction
23		Human SARS-coronavirus Nucleoprotein	25 ng/ml	Cross reaction
24		MERS-CoV Nucleoprotein	0.25 mg/ml	No cross reaction
25		Human Metapneumovirus (hMPV) 16 Type A1	$1.06 \times 10^{6.0}$ PFU/ml	No cross reaction
26		Adenovirus C1 Ad. 71	1.54×10^7 PFU/ml	No cross reaction
27		Enterovirus (EV68)	1.97×10^7 PFU/ml	No cross reaction

No.	Types of Specimen	Cross Reaction Substance	Final Test Concentration	Test Result
1	Other Microorganism	<i>Staphylococcus aureus</i>	$1.0 \times 10^{6.0}$ CFU/ml	No cross reaction
2		<i>Staphylococcus saprophyticus</i>	$1.0 \times 10^{6.0}$ CFU/ml	No cross reaction
3		<i>Neisseria</i> sp. (<i>Neisseria lactamica</i>)	$1.0 \times 10^{6.0}$ CFU/ml	No cross reaction

ENGLISH

No.	Types of Specimen	Cross Reaction Substance	Final Test Concentration	Test Result
4	Other Microorganism	<i>Escherichia coli</i>	$1.0 \times 10^{6.0}$ CFU/ml	No cross reaction
5		<i>Staphylococcus haemolyticus</i>	$1.0 \times 10^{6.0}$ CFU/ml	No cross reaction
6		<i>Streptococcus pyogenes</i>	$1.0 \times 10^{6.0}$ CFU/ml	No cross reaction
7		<i>Streptococcus salivarius</i>	$1.0 \times 10^{6.0}$ CFU/ml	No cross reaction
8		<i>Hemophilus parahaemolyticus</i>	$1.0 \times 10^{6.0}$ CFU/ml	No cross reaction
9		<i>Proteus vulgaris</i>	$1.0 \times 10^{6.0}$ CFU/ml	No cross reaction
10		<i>Moraxella catarrhalis</i>	$1.0 \times 10^{6.0}$ CFU/ml	No cross reaction
11		<i>Klebsiella pneumoniae</i>	$1.0 \times 10^{6.0}$ CFU/ml	No cross reaction
12		<i>Fusobacterium necrophorum</i>	$1.0 \times 10^{6.0}$ CFU/ml	No cross reaction
13		<i>Mycobacterium tuberculosis</i>	10 mg/ml	No cross reaction
14		Pooled human nasal wash	N/A*	No cross reaction
15		<i>Chlamydia pneumoniae</i>	9.1×10^7 CFU/ml	No cross reaction
16		<i>Haemophilus influenzae</i>	3.4×10^8 CFU/ml	No cross reaction
17		<i>Legionella pneumophila</i>	1.2×10^6 CFU/ml	No cross reaction
18		<i>Streptococcus pneumoniae</i> (4752-98 [Maryland (D1)6B-17])	$\geq 10^6$ CFU/ml	No cross reaction
19		<i>Bordetella pertussis</i>	4.4×10^9 CFU/ml	No cross reaction

ENGLISH

No.	Types of Specimen	Cross Reaction Substance	Final Test Concentration	Test Result
20	Other Microorganism	<i>Mycoplasma pneumoniae</i> (Mutant 22)	1.6 X 10 ⁹ CFU/ml	No cross reaction
21		<i>Pneumocystis jirovecii</i> (PJP)	1.0 X 10 ⁸ nuclei/ml	No cross reaction

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No.	Types of Specimen	Interfering Substances	Final Test Concentration	Test Result
1	Endogenous Substance	Mucin	0.5%	No Interference
2		Hemoglobin	100 mg/L	No Interference
3		Triglycerides	1.5 mg/L	No Interference
4		Icteric (Bilirubin)	40 mg/dL	No Interference
5		Rheumatoid factor	200 IU/ml	No Interference
6		Anti-nuclear antibody	>1:40	No Interference
7		Pregnant	10-fold dilution	No Interference
8	Exogenous Substance	Guaiacol glyceryl ether	1 µg/ml	No Interference
9		Albuterol	0.005 mg/dL	No Interference
10		Ephedrine	0.1 mg/ml	No Interference
11		Chlorpheniramine	0.08 mg/dL	No Interference
12		Diphenhydramine	0.08 mg/dL	No Interference
13		Ribavirin	26.7 µg /ml	No Interference
14		Oseltamivir	0.04 mg/dL	No Interference
15		Zanamivir	17.3 µg /ml	No Interference
16		Phenylephrine hydrochloride	15% v/v	No Interference
17		Oxymetazolin hydrochloride	15% v/v	No Interference
18		Amoxicillin	5.4 mg/dL	No Interference

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
ENGLISH

PREPARATION

- 1 Allow all kit components to reach a temperature between 15-30°C prior to testing for 30 minutes.
Note: Healthcare professionals should comply with personal safety guidelines including the use of personal protective equipment.
- 2 **Open the package and look for the following:**
 1. Test device with desiccant in individual foil pouch
 2. Buffer
 3. Extraction tube
 4. Extraction tube cap
 5. Positive control swab
 6. Negative control swab
 7. Sterilized nasopharyngeal swabs for sample collection
 8. Tube rack
 9. Quick reference guide (Nasopharyngeal)
 10. Instructions for use
- 3 Carefully read these instructions prior to using Panbio™ COVID-19 Ag Rapid Test Device kit.
- 4 Look at the expiration date of the kit box. If the expiration date has passed, use another kit.
- 5 **Open the foil pouch and look for the following:**
 1. Result window
 2. Specimen well

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


- 3 Tilt the patient's head back. Insert the swab through the nostril. Gently rub and roll the swab, 3-4 times. Leave the swab in place for several seconds. Slowly remove swab.

- 4 Insert the swab specimen in the extraction tube. Swirl the swab tip in the buffer fluid inside the extraction tube, pushing into the wall of the extraction tube at least five times and then squeeze out the swab by squeezing the extraction tube with your fingers.

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


⚠ **Caution:** The presence of any test line (T), no matter how faint, indicates a positive result.

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

1. Rothen HA, Byrareddy SN. The epidemiology and pathogenesis of coronavirus disease (COVID-19) outbreak. *J Autoimmun.* 2020; Feb 26:102433. doi:10.1016/j.jaut.2020.102433.
2. Guo YR, Cao QD, Hong ZS, et al. The origin, transmission and clinical therapies on coronavirus disease 2019 (COVID-19) outbreak-an update on the status. *Mil Med Res.* 2020; Mar 13; 7(1):11. doi:10.1186/s40779-020-00240-0.
3. Lai CC, Shih TP, Ko WC, et al. Severe acute respiratory syndrome coronavirus 2 (SARS-Cov-2) and coronavirus disease-2019 (COVID-19): The epidemic and the challenges. *Int J Antimicrob Agents.* 2020; Mar 55(3): 105924. doi: 10.1016/j.ijantimicag.2020.105924.
4. *In Vitro Diagnostic Assays for COVID-19: Recent Advances and Emerging Trends* (Sandeep Kumar Vashist, 2020 April 05: diagnostics)
5. Nano Research for COVID-19 (<https://dx.doi.org/10.1021/acsnano.0c02540>)
6. Coronavirus (COVID-19) Update: FDA Authorizes First Antigen Test to Help in the Rapid Detection of the Virus that Causes COVID-19 in Patients (Stephen M. Hahn M.D. 2020 May 09: Commissioner of Food and Drugs)
7. Current Intelligence Bulletin 13: Explosive Aside Hazard DHHS (NIOSH) Publication Number 78-127 August 16, 1976
8. CDC. Discontinuation of Transmission-Based Precautions and Disposition of Patients with COVID-19 in Healthcare Settings (Interim Guidance). (2020).
9. CDC. Duration of Isolation and Precautions for Adults with COVID-19. (2020).
10. Bullard, et al. Predicting Infectious Severe Acute Respiratory Syndrome Coronavirus 2 From Diagnostic Samples. *CID.* 2020; Nov 15;71 (10); DOI: 10.1093/cid/cia638

**GLOSSARY OF SYMBOLS /
SYMBOLVERZEICHNIS / GLOSARIO DE
SÍMBOLOS / GLOSSAIRE DES SYMBOLES /
GLOSSARIO DEI SIMBOLI / GLOSSÁRIO DE
SÍMBOLOS / ГЛОССАРИЙ ОБОЗНАЧЕНИЙ**

	Temperature limitation Temperaturbegrenzung Limitación de temperatura Limitation de température Limitazione di temperatura Limitação de temperatura Температурный диапазон
	For in vitro diagnostic use only Medizinprodukt für in-vitro Diagnostik Sólo para uso diagnóstico in vitro Pour un usage de diagnostic in vitro uniquement Ad uso esclusivo diagnostico in vitro Somente para uso para diagnóstico in vitro Медицинское изделие для диагностики In Vitro
	Do not reuse Nicht wiederverwenden No reutilizar Ne pas réutiliser Non riutilizzare Não reutilizar Запрет на повторное использование
	Do not use if package is damaged Bei beschädigter Verpackung nicht verwenden No lo use si el paquete está dañado Ne pas utiliser si le colis est endommagé Non utilizzare se la confezione è danneggiata Não use se o pacote estiver danificado Не использовать при повреждении упаковки

	Lot Number Chargencode Número de lote Numéro de lot Numero di lotto Número de Lote Номер серии
	Catalog Number Artikelnummer Número de catalogo Numéro de catalogue Numero di catalogo Número no Catálogo Номер изделия по каталогу
	Consult instructions for use Gebrauchsanleitung beachten Consultar instrucciones de uso Consulter les instructions d'utilisation Consultare le istruzioni per l'uso Consulte as instruções de uso Обратитесь к инструкции по применению
	Keep dry Trocken aufbewahren Mantener seco Garder au sec Mantenere asciutto Manter seco Береечь от влаги
	Biological Risks Biologisches Risiko Riesgos biológicos Risques biologiques Rischi biologici Riscos Biológicos Биологическая опасность
	Use By Verwendbar bis Usar por Utiliser par Utilizzare per Usar até Использовать до

	Manufacturer Hersteller Fabricante Fabricant Produttore Fabricante Производитель
	Date of manufacture Herstellungsdatum Fecha de manufactura Date de fabrication Data di produzione Data de fabricaço Дата производства
	Keep away from sunlight Von Sonnenlicht fernhalten Mantener alejado de la luz solar Tenir à l'écart de la lumière du soleil Tenere lontano dalla luce solare Manter longe da luz solar Не допускать воздействия солнечного света
	CE mark CE Zeichen Marca CE Marquage CE Marcatura CE Marca CE Знак соответствия продукции техническим регламентам ЕС
	Contains sufficient for X tests Ausreichend für X Prüfungen Contiene suficiente para X pruebas Contient suffisamment pour les tests X Contenuto sufficiente per X test Contém suficiente para X testes Содержимого достаточно для проведения n-количества тестов
	Caution Achtung Precaución Attention Attenzione Cuidado Осторожно! Обратитесь к инструкции по применению

STERILE EO	Sterilized using ethylene oxide Sterilisiert mit Ethylenoxid Esterilizado con óxido de etileno. Stérilisé à l'oxyde d'éthylène Sterilizzato con ossido di etilene Esterilizado com óxido de etileno Стерилизовано оксидом этилена
STERILE R	Sterilized using irradiation Sterilisiert durch Bestrahlung Esterilizado mediante irradiación. Stérilisé par irradiation Sterilizzato con irradiazione Esterilizado por irradiação Радиационная стерилизация
	Do not re-sterilize Nicht sterilisieren No volver a esterilizar Ne pas re-stériliser Non risterilizzare Não reesterilize Не стерилизовать повторно
CONTROL -	Negative control Negative Kontrolle Control negativo Contrôle négatif Controllo negativo Controle negativo Отрицательный контроль
CONTROL +	Positive control Positive Kontrolle Control positivo Contrôle positif Controllo positivo Controle positivo Положительный контроль

TECHNICAL SUPPORT / TECHNISCHER SUPPORT / SOPORTE TÉCNICO / SUPPORT TECHNIQUE / SUPPORTO TECNICO / SUPORTE TÉCNICO / ОТДЕЛ ТЕХНИЧЕСКОЙ ПОДДЕРЖКИ	
Europe & Middle East / Europa & Mittlerer Osten / Europa y Medio Oriente / Europe & Moyen Orient / Europe & Middle East / Europa e Oriente Médio / Европа, Ближний Восток	+44 161 483 9032 EME.techsupport@abbott.com
Africa / Afrika / África / Afrique / Africa / África / Африка	+27 10 500 9700 ARCIS.techsupport@abbott.com
Russia & CIS / Russland & CIS / Rusia y CEI / Russie & CIS / Russia & CIS / Rússia e CEI / Россия и СНГ	+44 161 483 9032 ARCIS.techsupport@abbott.com
Asia Pacific / Asien Pazifik / Asia Pacífico / Asie Pacifique / Asia Pacific / Ásia-Pacífico / Азия, Океания	+61 7 3363 7100 AP.techsupport@abbott.com
Latin America / Lateinamerika / America Latina / Amérique Latine / Latin America / América Latina / Латинская Америка	+57 1 482 4033 LA.techsupport@abbott.com
Canada/Kanada/Canada/ Canada/Canadá/Canadá/ Канада	+1 800 818 8335 Canproductsupport@abbott.com

- MEMO -

- 142 -

 **Abbott Rapid Diagnostics Jena GmbH**
Ortenberg 1, D-07743 Jena, Germany
abbott.com/post

Date Issued : 2021.03
41FK10/41FK20-07-A2

© 2021 Abbott. All rights reserved.
All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

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
TECHNICAL SUPPORT: ABBOTT.COM/POCT

© 2020 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

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
Orlweg 1 D-07743 Jena, Germany
Tel. +49(0)3641-31110 Fax +49(0)3641-311120


File-No.: 1637 / 2022 / L

I hereby notarise that the document bears Mr. Singer signature. This was signed in my presence by

Mr. Peter Singer,
born on 19th May 1987,
business address: Orlaweg 1, 07743 Jena
who identified himself by means of his
identity papers.

Mr. Singer declares not to act on his own behalf, but as Regulatory Affairs Manager for the German company "Abbott Rapid Diagnostics Jena GmbH" with its registered office in Jena, registered at the Court of Jena, under Commercial Register No. HRB 207136.

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Signed: Peter Singer
Peter Singer,
Regulatory Affairs Specialist,
Abbott Rapid Diagnostics Jena GmbH

Date: 16.12.2022

Abbott Rapid Diagnostics Jena GmbH

Orlaweg 1 D-07743 Jena, Germany
Tel. +49(0)3641-31110 Fax +49(0)3641-3111120

Abbott Rapid Diagnostics Jena GmbH
Orlaweg 1
D-07743 Jena / Germany
Phone: +49 (0) 3641 3111-0
Fax.: +49 (0) 3641 3111-120
E-Mail: atj.info@abbott.com
Web: www.abbott.com/poct

Managing Directors/ Geschäftsführer:
Dr. Hans Olaf Such
John Arthur McCoy Jr.
Benjamin Edward Oosterbaan
Damian Patrick Halloran

Company reg.: HRB 207136 AG Jena
VAT No.: DE194889265

Commerzbank AG, Jena
Bank code: 82040000, Acc.: 258882000
IBAN: DE08 8204 0000 0258 8820 00
SWIFT (BIC): COBADEFF821
Sparkasse Jena
Bank code: 83053030, Acc.: 142093
IBAN: DE83 8305 3030 0000 1420 93
SWIFT (BIC): HELADEF1JEN

REF 41FK11/41FK21

Panbio-
**COVID-19 Ag Rapid
Test Device**
(NASAL)

In vitro diagnostic rapid test for qualitative detection of SARS-CoV-2 antigen (Ag)

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EN

Introduction

The Coronavirus disease (COVID-19) is an infectious disease caused by a newly discovered coronavirus, severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2). The SARS-CoV-2 is a β -coronavirus, which is an enveloped non-segmented positive-sense RNA virus¹. It is spread by human-to-human transmission via droplets or direct contact, and infection has been estimated to have a mean incubation period of 6.4 days and a basic reproduction number of 2.24-3.58. Among patients with pneumonia caused by SARS-CoV-2, fever was the most common symptom, followed by cough². The main IVD assays used for COVID-19 employ real-time reverse transcriptase-polymerase chain reaction (RT-PCR) that takes a few hours³. The availability of a cost-effective, rapid point-of-care diagnostic test is critical to enable healthcare professionals to aid in the diagnosis of patients and prevent further spread of the virus³. Antigen tests will play a critical role in the fight against COVID-19⁴.

Test Principle

Panbio™ COVID-19 Ag Rapid Test Device contains a membrane strip, which is pre-coated with immobilized anti-SARS-CoV-2 antibody on the test line and mouse monoclonal anti-chicken IgY on the control line. Two types of conjugates (human IgG specific to SARS-CoV-2 Ag gold conjugate (binds to the nucleocapsid protein) and chicken IgY gold conjugate) move upward on the membrane chromatographically and react with anti-SARS-CoV-2 antibody and pre-coated mouse monoclonal anti-chicken IgY respectively. For a positive result, human IgG specific to SARS-CoV-2 Ag gold conjugate and anti-SARS-CoV-2 antibody will form a test line in the result window. Neither the test line nor the control line are visible in the result window prior to applying the patient specimen. A visible control line is required to indicate a test result is valid.

Intended Use

Panbio™ COVID-19 Ag Rapid Test Device is an *in vitro* diagnostic rapid test for the qualitative detection of SARS-CoV-2 antigen (Ag) in human nasal swab specimens from individuals who meet COVID-19 clinical and / or epidemiological criteria. Panbio™ COVID-19 Ag Rapid Test Device is for professional use only and is intended to be used as an aid in the diagnosis of SARS-CoV-2 infection. The product may be used in any laboratory and non-laboratory environment that meets the requirements specified in the Instructions for Use and local regulation. The test provides preliminary test results. Negative results don't preclude SARS-CoV-2 infection and they cannot be used as the sole basis for treatment or other management decisions. Negative results must be combined with clinical observations, patient history, and epidemiological information. The test is not intended to be used as a donor screening test for SARS-CoV-2.

ENGLISH

Kit Variants

- 41FK11 No 2D barcode printed on the contained test devices
- 41FK21 Contains test devices with a 2D barcode printed on the test device, which encodes traceability information for the product

Materials Provided

- 25 Test devices with desiccant in individual foil pouch
- Buffer (1 x 9 ml/bottle)
- 25 Extraction tubes
- 25 Extraction tube caps
- 1 Positive control swab
- 1 Negative control swab
- 25 Sterilized nasal swabs for sample collection
- 1 Tube rack
- 1 Quick Reference Guide (Nasal)
- 1 Instructions for use

Materials Required but not Provided

- Personal Protective Equipment per local recommendations (i.e. gown/lab coat, face mask, face shield/eye goggles and gloves), Timer, Biohazard container

Active Ingredients of Main Components

- 1 Test device Gold conjugate: Human IgG specific to SARS-CoV-2 Ag gold colloid and Chicken IgY - gold colloid, Test line: Mouse monoclonal anti-SARS-CoV-2, Control line: Mouse monoclonal anti-Chicken IgY
- Buffer Tricine, Sodium Chloride, Tween 20, Sodium Azide (<0.1%), Proclin 300

Storage and Stability

1. The test kit should be stored at a temperature between 2-30 °C. Do not freeze the kit or its components.
Note: When stored in a refrigerator, all kit components must be brought to room temperature (15-30 °C) for a minimum of 30 minutes prior to performing the test. Do not open the pouch while components come to room temperature.
2. The Buffer bottle may be opened and resealed for each assay. The Buffer cap should be firmly sealed between each use. The Buffer is stable until expiration date if kept at 2-30 °C.
3. Perform the test immediately after removing the test device from the foil pouch.
4. Do not use the test kit beyond its expiration date.

ENGLISH

5. The shelf life of the kit is as indicated on the outer package.
6. Do not use the test kit if the pouch is damaged or the seal is broken.
7. Direct swab specimens should be tested immediately after collection. If immediate testing is not possible, the swab specimen can be kept in an extraction tube filled with extraction buffer (300 µl) at room temperature (15-30 °C) for up to two hours prior to testing.

Warnings

1. For in vitro diagnostic use only. Do not reuse the test device and kit components.
2. These instructions must be strictly followed by a trained healthcare professional to achieve accurate results. All users have to read the instruction prior to performing a test.
3. Do not eat or smoke while handling specimens.
4. Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
5. Avoid splashing or aerosol formation of specimen and buffer.
6. Clean up spills thoroughly using an appropriate disinfectant.
7. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials (i.e. swab, extraction tube, test device) in a biohazard container as if they were infectious waste and dispose according to applicable local regulations.
8. Do not mix or interchange different specimens.
9. Do not mix reagent of different lots or those for other products.
10. Do not store the test kit in direct sunlight.
11. To avoid contamination, do not touch the head of provided swab when opening the swab pouch.
12. The sterilized swabs should be used only for nasal specimen collection.
13. To avoid cross-contamination, do not reuse the sterilized swabs for specimen collection.
14. Do not dilute the collected swab with any solution except for the provided extraction buffer.
15. The buffer contains <0.1% sodium azide as a preservative which may be toxic if ingested. When disposed of through a sink, flush with a large volume of water.⁷
16. Do not use the positive or negative control swab for specimen collection.

Test Procedure (Refer to Figure)

Nasal swab Specimens

Note: Healthcare professionals should comply with personal safety guidelines including the use of personal protective equipment.

ENGLISH

Test Preparation

1. Allow all kit components to reach a temperature between 15-30 °C prior to testing for 30 minutes.
2. Remove the test device from the foil pouch prior to use. Place on a flat, horizontal and clean surface.
3. Hold the buffer bottle vertically and fill the extraction tube with buffer fluid until it flows up to the Fill-line of the extraction tube (300 µl).
⚠ Caution: If the amount of buffer is excessive or insufficient, an improper test result may occur.
4. Place the extraction tube in the tube rack.

Nasal Mid-Turbinate (NMT) Specimen Collection & Extraction

Specimens are collected by the professional user as described below. Alternatively, nasal specimen collection steps 1-3 can be completed by the patient according to oral instructions and under supervision of the professional user. For supervised patient self-collection, the swab is handed to the patient by the professional user and after sampling, the patient hands the swab back to the professional user to complete the remaining steps of the procedure.

1. Tilt the patient's head back 70 degrees. While gently rotating the swab, insert swab less than one inch (about 2 cm) into nostril (until resistance is met at the turbinates).
⚠ Caution: Ensure that the patient's head is kept still during nasal specimen collection, as sudden movements may cause swab stick breakage.
2. Rotate the swab five times against the nasal wall then slowly remove from the nostril.
3. Using the same swab repeat the collection procedure with the second nostril.
Note: Ensure a minimum waiting period of 24 hours before performing a new nasal sampling from both nostrils (e.g. for a repeat test).
⚠ Caution: If the nasal swab stick breaks prior to obtaining a nasal specimen, repeat specimen collection with a new swab. If the nasal swab breaks during sampling, consultation with a medical healthcare professional is recommended to determine and initiate necessary treatment and monitoring.
4. Swirl the swab tip in the buffer fluid inside the extraction tube, pushing into the wall of the extraction tube at least five times and then squeeze out the swab by squeezing the extraction tube with your fingers.
5. Break the swab at the breakpoint and close the cap of extraction tube.

Reaction with Test Device

1. Open the dropping nozzle cap at the bottom of the extraction tube.
2. Dispense 5 drops of extracted specimens vertically into the specimen well (S) on the device. Do not handle or move the test device until the test is complete and ready for reading.
⚠ Caution: Bubbles that occur in the extraction tube can lead to inaccurate

ENGLISH

results. If you are unable to create sufficient drops, this may be caused by clogging in the dispensing nozzle. Shake the tube gently to release the blockage until you observe free drop formation.

3. Close the nozzle and dispose of the extraction tube containing the used swab according to your local regulations and biohazard waste disposal protocol.
4. Start timer. Read result at 15 minutes. Do not read results after 20 minutes.
5. Dispose of the used device according to your local regulations and biohazard waste disposal protocol.

Positive / Negative Control Swab

⚠ Caution: Control use only. Do not use the positive or negative control swab for specimen collection.

Note: Please refer to the External Quality Control section of this Instructions for use for the frequency of testing external quality control swabs.

1. Hold the buffer bottle vertically and fill the extraction tube with buffer fluid until it flows up to the Fill-line of the extraction tube (300 µl).
⚠ Caution: If the amount of buffer is excessive or insufficient, an improper test result may occur.
2. Place the extraction tube in the tube rack.
3. Insert the positive or negative control swab in the buffer fluid inside of the extraction tube and soak the swab for 1 minute. Swirl the control swab tip in the buffer fluid inside of the extraction tube, pushing into the wall of the extraction tube at least five times and then squeeze out the swab by squeezing the extraction tube with your fingers.
4. Dispose of the used control swab in accordance with your biohazard waste disposal protocol.
5. Close the cap of the extraction tube.
6. Follow the above test procedure [Reaction with Test Device].

Test Interpretation (Refer to Figure)

1. **Negative result:** The presence of only the control line (C) and no test line (T) within the result window indicates a negative result.
2. **Positive result:** The presence of the test line (T) and the control line (C) within the result window, regardless of which line appears first, indicates a positive result.
⚠ Caution: The presence of any test line (T), no matter how faint, indicates a positive result.
3. **Invalid result:** If the control line (C) is not visible within the result window after performing the test, the result is considered invalid.

ENGLISH

Test Limitations

1. The contents of this kit are to be used for the professional and qualitative detection of SARS-CoV-2 antigen from nasal swab. Other specimen types may lead to incorrect results and must not be used.
2. Failure to follow the instructions for test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
3. A negative test result may occur if the specimen was collected, extracted or transported improperly. A negative test result does not eliminate the possibility of SARS-CoV-2 infection and should be confirmed by viral culture or a molecular assay.
4. Positive test results do not rule out co-infections with other pathogens.
5. Test results must be evaluated in conjunction with other clinical data available to the physician.
6. Reading the test results earlier than 15 minutes or later than 20 minutes may give incorrect results.
7. Panbio™ COVID-19 Ag Rapid Test Device is not intended to detect from defective (non-infectious) virus during the later stages of viral shedding that might be detected by PCR molecular tests.⁸
8. Positive results may occur in cases of infection with SARS-CoV.

Quality Control

1. Internal Quality Control:

The test device has a test line (T) and a control line (C) on the surface of the test device. Neither the test line nor the control line are visible in the result window before applying a specimen. The control line is used for procedural control and should always appear if the test procedure is performed properly and the test reagents of the control line are working.

2. External Quality Control:

The controls are specifically formulated and manufactured to ensure performance of the Panbio™ COVID-19 Ag Rapid Test Device and are used to verify the user's ability to properly perform the test and interpret the results. The Positive Control contains recombinant SARS-CoV-2 nucleocapsid protein, which is not contagious. The Positive Control will produce a positive test result and has been manufactured to produce a visible test line (T). The Negative Control will produce a negative test result. Control swabs are not specific for a particular Panbio™ COVID-19 Ag Rapid Test Device lot and may be used between test device lots until the swabs' expiry dates.

Good laboratory practice suggests the use of positive and negative controls to ensure that:

- Test reagents are working, and
- The test is correctly performed.

The external controls can be run under any of the following circumstances:

- By a new operator prior to performing testing on patient specimens,
- When receiving a new test shipment,
- At periodic intervals as dictated by local requirements, and/or by the user's Quality Control procedures.

Performance Characteristics

1. External evaluation of Panbio™ COVID-19 Ag Rapid Test Device (Symptomatic)

Clinical performance of Panbio™ COVID-19 Ag Rapid Test Device was determined by testing 104 positive nasal swab specimens and 404 negative specimens for SARS-CoV-2 antigen (Ag) to have a sensitivity of 98.1% (95% CI: 93.2-99.8%) and a specificity of 99.8% (95% CI: 98.6-100.0%). Clinical specimens were determined to be positive or negative using an FDA EUA RT-PCR reference method. The individuals on which the reported sensitivity and specificity are based also had a nasopharyngeal swab taken, which was tested in the FDA EUA approved RT-PCR.

Panbio™ COVID-19 Ag Rapid Test Device Results

		Nasal PCR Test Result		
		Positive	Negative	Total
Panbio™ COVID-19 Ag Rapid Test Device Result (nasal swab specimens)	Positive	102	1	103
	Negative	2	403	405
	Total	104	404	508
		Sensitivity	Specificity	Overall Percent Agreement
		98.1% [93.2%; 99.8%]	99.8% [98.6%; 100.0%]	99.4% [98.3%; 99.9%]

- Performance data was calculated from a study of individuals suspected of exposure to COVID-19 or who have presented with symptoms in the last 7 days.
- Stratification of the positive specimens post onset of symptoms or suspected exposure between 0-3 days has a sensitivity of 100.0% (95% CI: 92.3-100.0%; n=46) and 4-7 days has a sensitivity of 96.6% (95% CI: 88.1-99.6%; n=58).
- Positive agreement of the Panbio™ COVID-19 Ag Rapid Test Device is higher with samples of Ct values ≤30 with a sensitivity of 100.0% (95% CI: 96.0-100.0%) and Ct values ≤33 with a sensitivity of 99.0% (95% CI: 94.5-100.0%). As indicated in References 8-10, patients with Ct value >30 are no longer contagious.^{8, 9, 10}
- The clinical performance data was also calculated vs nasopharyngeal swab specimens using an FDA EUA RT-PCR reference and has a sensitivity of 91.1% (95% CI: 84.2-95.6%) and specificity of 99.7% (95% CI: 98.6-100.0%).

2. External evaluation of Panbio™ COVID-19 Ag Rapid Test Device (Asymptomatic)

Clinical performance of Panbio™ COVID-19 Ag Rapid Test Device was determined by testing 483 asymptomatic subjects for SARS-CoV-2 antigen (Ag). Clinical specimens were determined to be positive or negative using an FDA EUA RT-PCR reference method.

The positive results (n=50) were stratified by the comparator method cycle threshold (Ct) counts and assessed to better understand the correlation of product performance, as a surrogate for the amount of virus present in the clinical sample. A lower Ct value corresponds to a higher virus concentration. As presented in the table below, the positive agreement increases with lower Ct values.

The specificity (n=433) was 100% with 95% CI [99.2%; 100.0%].

The results for sensitivity are summarized in the following table:

	All Nasal PCR Positive Samples (n=50)	Ct values ≤ 33 (n=40)	Ct values ≤ 30 Ct (n=32)
Sensitivity [CI 95%]	66.0% [51.2%; 78.8%]	80.0% [64.4%; 90.9%]	93.8% [79.2%; 99.2%]

As indicated in References 8-10, patients with Ct value >30 are no longer contagious.^{8,9,10}

3. External evaluation of Panbio™ COVID-19 Ag Rapid Test Device (Self-Collected Swab)

The clinical performance of Panbio™ COVID-19 Ag Rapid Test Device was assessed in 287 symptomatic subjects (≥16 years of age) who collected their swab specimen (self swabbing) under the direction and supervision of a trained professional. The swab was then handed to the trained professional who executed the remaining steps of the procedure. The trained professional also collected a nasopharyngeal swab from each subject to be used as a reference specimen. The reference specimen was tested on the Panbio™ COVID-19 Ag Rapid Test Device.

The results are summarized in the following table:

		Panbio™ COVID-19 Ag Rapid Test Device (Nasopharyngeal)		
		Positive	Negative	Total
Panbio™ COVID-19 Ag Rapid Test Device (Nasal) – Self-collected Swab	Positive	110	0	110
	Negative	2	175	177
	Total	112	175	287
		Positive Agreement	Negative Agreement	Overall Percent Agreement
		98.2% [93.7%; 99.8%]	100.0% [97.9%; 100.0%]	99.3% [97.5%; 99.9%]

4. Detection Limit

Panbio™ COVID-19 Ag Rapid Test Device was confirmed to detect $2.5 \times 10^{1.8}$ TCID₅₀/ml of SARS-CoV-2 which was isolated from a COVID-19 confirmed patient in Korea.

5. Hook Effect

There is no hook effect at $1.0 \times 10^{5.8}$ TCID₅₀/ml of SARS-CoV-2 which was isolated from a COVID-19 confirmed patient in Korea.

6. Cross Reactivity

Cross-reactivity of Panbio™ COVID-19 Ag Rapid Test Device was evaluated by testing 46 viruses and 21 other microorganisms. The final test concentrations of viruses and other microorganisms are documented in the Table below. The following viruses and other microorganisms except the Human SARS-coronavirus Nucleoprotein have no effect on the test results of Panbio™ COVID-19 Ag Rapid Test Device.

Panbio™ COVID-19 Ag Rapid Test Device has cross-reactivity with Human SARS-coronavirus Nucleoprotein at a concentration of 25 ng/ml or more because SARS-CoV has high homology (79.6%) to the SARS-CoV-2.

No.	Types of Specimen	Cross Reaction Substance	Final Test Concentration	Test Result
1	Virus	Adenovirus Type 1	1.54×10^7 PFU/ml	No cross reaction
2		Adenovirus Type 5	4.0×10^8 PFU/ml	No cross reaction
3		Adenovirus Type 7	2.0×10^9 PFU/ml	No cross reaction
4		Enterovirus (EV68)	2.0×10^7 PFU/ml	No cross reaction
5		Echovirus2	$7.0 \times 10^{5.5}$ PFU/ml	No cross reaction

ENGLISH

No.	Types of Specimen	Cross Reaction Substance	Final Test Concentration	Test Result
6	Virus	Echovirus11	$3.5 \times 10^{6.25}$ PFU/ml	No cross reaction
7		Enterovirus D68	2.0×10^7 PFU/ml	No cross reaction
8		Human herpesvirus (HSV) 1	$3.5 \times 10^{7.5}$ PFU/ml	No cross reaction
9		Human herpesvirus (HSV) 2	$3.5 \times 10^{5.75}$ PFU/ml	No cross reaction
10		Mumps Virus Ag	1.1×10^5 PFU/ml	No cross reaction
11		Influenza virus A (H1N1) Strain (A/Virginia/ATCC1/2009)	2.6×10^5 PFU/ml	No cross reaction
12		Influenza virus A (H1N1) Strain (A/WS/33)	$3.5 \times 10^{7.25}$ PFU/ml	No cross reaction
13		Influenza virus A(H1N1) Strain (A/California/08/2009/pdm09)	1.1×10^8 PFU/ml	No cross reaction
14		Influenza virus B Strain (B/Lee/40)	$3.5 \times 10^{6.25}$ PFU/ml	No cross reaction
15		Parainfluenza Type 1	2.1×10^8 PFU/ml	No cross reaction
16		Parainfluenza Type 2	3.5×10^5 PFU/ml	No cross reaction
17		Parainfluenza Type 3	4.6×10^7 PFU/ml	No cross reaction
18		Parainfluenza Type 4A	2.0×10^7 PFU/ml	No cross reaction
19		Respiratory syncytial virus (RSV) type A	3.0×10^5 PFU/ml	No cross reaction
20		Respiratory syncytial virus (RSV) type B	3.9×10^5 PFU/ml	No cross reaction
21		Rhinovirus A16	8.8×10^5 PFU/ml	No cross reaction
22		HCoV-HKU1	1.5mg/ml	No cross reaction
23		HCoV-NL63	1.2×10^5 PFU/ml	No cross reaction
24		HCoV-OC43	6.2×10^5 PFU/ml	No cross reaction

No.	Types of Specimen	Cross Reaction Substance	Final Test Concentration	Test Result
25	Virus	HCoV-229E	1.1×10^6 PFU/ml	No cross reaction
26		Human SARS-coronavirus Nucleoprotein	25 ng/ml	Cross Reaction
27		MERS-CoV Nucleoprotein	0.25 mg/ml	No cross reaction
28		Human Metapneumovirus (hMPV) 16 Type A1	1.1×10^6 PFU/ml	No cross reaction
29		Adenovirus Type 2	1.96×10^7 PFU/ml	No cross reaction
30		Adenovirus Type 3	$1.4 \times 10^{6.5}$ PFU/ml	No cross reaction
31		Adenovirus Type 4	$3.5 \times 10^{6.5}$ PFU/ml	No cross reaction
32		Enterovirus C	6.0×10^7 PFU/ml	No cross reaction
33		Influenza virus A(H3N2) Strain (A/Hong Kong/8/68)	$3.5 \times 10^{5.5}$ PFU/ml	No cross reaction
34		Influenza virus A(H5N1)	1.5 mg/ml	No cross reaction
35		Influenza virus B Strain (Victoria)	5.46×10^6 PFU/ml	No cross reaction
36		Rhinovirus 14	1.6×10^8 PFU/ml	No cross reaction
37		Human cytomegalovirus	7.0×10^5 PFU/ml	No cross reaction
38		Norovirus	7.14×10^7 PFU/ml	No cross reaction
39		Varicella-zoster virus	1.96×10^4 PFU/ml	No cross reaction
40		Measles virus	6.1×10^5 PFU/ml	No cross reaction
41		EB virus	5.6×10^8 copies/ml	No cross reaction
42		Influenza virus(H7N9)	1.5mg/ml	No cross reaction

ENGLISH

No.	Types of Specimen	Cross Reaction Substance	Final Test Concentration	Test Result
43	Virus	Influenza virus B Strain (Yamagata)	2.73×10^{10} PFU/ml	No cross reaction
44		Rhinovirus 54	$3.5 \times 10^{5.67}$ PFU/ml	No cross reaction
45		Rotavirus	1.12×10^7 PFU/ml	No cross reaction
46		Adenovirus type 11	3.0×10^6 PFU/ml	No cross reaction

No.	Types of Specimen	Cross Reaction Substance	Final Test Concentration	Test Result
1	Other Microorganism	<i>Staphylococcus saprophyticus</i>	7.9×10^7 CFU/ml	No cross reaction
2		<i>Neisseria sp. (Neisseria lactamica)</i>	6.8×10^8 CFU/ml	No cross reaction
3		<i>Staphylococcus haemolyticus</i>	1.4×10^{10} CFU/ml	No cross reaction
4		<i>Streptococcus salivarius</i>	7.84×10^7 CFU/ml	No cross reaction
5		<i>Hemophilus parahaemolyticus</i>	8.8×10^8 CFU/ml	No cross reaction
6		<i>Proteus vulgaris</i>	2.9×10^7 CFU/ml	No cross reaction
7		<i>Moraxella catarrhalis</i>	1.9×10^8 CFU/ml	No cross reaction
8		<i>Klebsiella pneumoniae</i>	2.0×10^7 CFU/ml	No cross reaction
9		<i>Fusobacterium necrophorum</i>	7.0×10^8 CFU/ml	No cross reaction
10		<i>Mycobacterium tuberculosis</i>	10mg/ml	No cross reaction
11		Pooled human nasal wash	N/A*	No cross reaction
12		<i>Streptococcus pyogenes</i>	3.6×10^7 CFU/ml	No cross reaction
13		<i>Mycoplasma pneumoniae</i>	4.0×10^8 CFU/ml	No cross reaction

ENGLISH

No.	Types of Specimen	Cross Reaction Substance	Final Test Concentration	Test Result
14	Other Microorganism	<i>Staphylococcus aureus</i>	1.3×10^8 CFU/ml	No cross reaction
15		<i>Escherichia coli</i>	6.8×10^6 CFU/ml	No cross reaction
16		<i>Chlamydia pneumoniae</i>	9.1×10^7 IFU/ml	No cross reaction
17		<i>Haemophilus influenzae</i>	3.4×10^8 CFU/ml	No cross reaction
18		<i>Legionella pneumophila</i>	1.2×10^6 CFU/ml	No cross reaction
19		<i>Streptococcus pneumoniae</i>	1.3×10^6 CFU/ml	No cross reaction
20		<i>Bordetella pertussis</i>	4.4×10^9 CFU/ml	No cross reaction
21		<i>Pneumocystis jirovecii</i> (PJP)	1.0×10^8 nuclei/ml	No cross reaction

* No concentration provided by supplier. Undiluted stock solution was tested.

7. Interfering Substances

The following 43 potentially interfering substances have no impact on Panbio™ COVID-19 Ag Rapid Test Device. The final test concentrations of the interfering substances are documented in the Table below.

No.	Types of Specimen	Interfering Substances	Final Test Concentration	Test Result
1	Endogenous Substance	Mucin	0.5%	No Interference
2		Hemoglobin	100 mg/L	No Interference
3		Triglycerides	1.5 mg/L	No Interference
4		Icteric (Bilirubin)	40 mg/dL	No Interference
5		Rheumatoid factor	200 IU/ml	No Interference
6		Anti-nuclear antibody	>1:40	No Interference
7		Pregnant	10-fold dilution	No Interference
8	Exogenous Substance	Guaiacol glyceryl ether	1 µg/ml	No Interference
9		Albuterol	0.005 mg/dL	No Interference
10		Ephedrine	0.1 mg/ml	No Interference
11		Chlorpheniramine	0.08 mg/dL	No Interference
12		Diphenhydramine	0.08 mg/dL	No Interference
13		Ribavirin	26.7 µg /ml	No Interference
14		Oseltamivir	0.04 mg/dL	No Interference

ENGLISH

No.	Types of Specimen	Interfering Substances	Final Test Concentration	Test Result
15	Exogenous Substance	Zanamivir	17.3 µg /ml	No Interference
16		Phenylephrine hydrochloride	15% v/v	No Interference
17		Oxymetazolin hydrochloride	15% v/v	No Interference
18		Amoxicillin	5.4 mg/dL	No Interference
19		Acetylsalicylic acid	3 mg/dL	No Interference
20		Ibuprofen	21.9 mg/dL	No Interference
21		Chlorothiazide	2.7 mg/dL	No Interference
22		Indapamide	140 ng/ml	No Interference
23		Glimepiride (Sulfonylureas)	0.164 mg/dL	No Interference
24		Acarbose	0.03 mg/dL	No Interference
25		Ivermectin	4.4 mg/L	No Interference
26		Lopinavir	16.4 µg/L	No Interference
27		Ritonavir	16.4 µg/L	No Interference
28		Chloroquine phosphate	0.99 mg/L	No Interference
29		Sodium chloride with preservatives	4.44 mg/ml	No Interference
30		Beclomethasone	4.79 ng/ml	No Interference
31		Dexamethasone	0.6 µg/ml	No Interference
32		Flunisolide	0.61 µg/ml	No Interference
33		Triamcinolone	1.18 ng/ml	No Interference
34		Budesonide	2.76 ng/ml	No Interference
35		Mometasone	1.28 ng/ml	No Interference
36		Fluticasone	2.31 ng/ml	No Interference
37		Sulfur	9.23 µg/ml	No Interference
38		Benzocaine	0.13 mg/ml	No Interference
39		Menthol	0.15 mg/ml	No Interference
40		Mupirocin	10 µg/ml	No Interference
41		Tobramycin	24.03 µg/ml	No Interference
42		Biotin	1.2 µg/ml	No Interference
43		HAMA	63.0 ng/ml	No Interference

ENGLISH

8. Repeatability & Reproducibility


Repeatability & Reproducibility of Panbio™ COVID-19 Ag Rapid Test Device was established using in-house reference panels containing negative specimens and a range of positive specimens. There were no differences observed within-run, between-run, between-lots, between-sites, and between-days.

ENGLISH

PREPARATION

- 1 Allow all kit components to reach a temperature between 15-30°C prior to testing for 30 minutes.
Note: Healthcare professionals should comply with personal safety guidelines including the use of personal protective equipment.
- 2 **Open the package and look for the following:**
 1. Test device with desiccant in individual foil pouch
 2. Buffer
 3. Extraction tube
 4. Extraction tube cap
 5. Positive control swab
 6. Negative control swab
 7. Sterilized nasal swabs for sample collection
 8. Tube rack
 9. Quick reference guide (Nasal)
 10. Instructions for use
- 3 Carefully read these instructions prior to using Panbio™ COVID-19 Ag Rapid Test Device kit.
- 4 Look at the expiration date of the kit box. If the expiration date has passed, use another kit.
- 5 **Open the foil pouch and look for the following:**
 1. Result window
 2. Specimen well

Then, label the device with the patient identifier.


••••• : 5 drops of the extracted specimen

TEST PROCEDURE

- 1 Hold the buffer bottle vertically and fill the extraction tube with buffer fluid until it flows up to the Fill-line of the extraction tube (300µl).

⚠ Caution: If the amount of buffer is excessive or insufficient, an improper test result may occur.


- 2 Place the extraction tube in the tube rack.

- 3 Tilt the patient's head back 70 degrees. While gently rotating the swab, insert swab less than one inch (about 2 cm) into nostril (until resistance is met at the turbinates). Rotate the swab five times against the nasal wall. Using the same swab repeat the collection procedure with the second nostril. Slowly remove swab from the nostril.


⚠ Caution: If the swab stick breaks during specimen collection, repeat specimen collection with a new swab.


ENGLISH

TEST PROCEDURE

- 4** Insert the swab specimen in the extraction tube. Swirl the swab tip in the buffer fluid inside the extraction tube, pushing into the wall of the extraction tube at least five times and then squeeze out the swab by squeezing the extraction tube with your fingers.

- 5** Break the swab at the breakpoint and close the cap of extraction tube.


- 6** Open the dropping nozzle cap at the bottom of the extraction tube.

- 7** Dispense 5 drops of extracted specimens vertically into the specimen well (S) on the device. Do not handle or move the test device until the test is complete and ready for reading.


⚠ Caution: Bubbles that occur in the extraction tube can lead to inaccurate results. If you are unable to create sufficient drops, this may be caused by clogging in the dispensing nozzle. Shake the tube gently to release the blockage until you observe free drop formation.

ENGLISH

TEST PROCEDURE

- 8** Close the nozzle and dispose of the extraction tube containing the used swab according to your local regulations and biohazard waste disposal protocol.

- 9** Start timer. Read result at 15 minutes. Do not read results after 20 minutes.

- 10** Dispose of the used device according to your local regulations and biohazard waste disposal protocol.


ENGLISH

TEST INTERPRETATION

NEGATIVE


The presence of only the control line (C) and no test line (T) within the result window indicates a negative result.

POSITIVE

The presence of the test line (T) and the control line (C) within the result window, regardless of which line appears first, indicates a positive result.

 **Caution:** The presence of any test line (T), no matter how faint, indicates a positive result.

INVALID

If the control line (C) is not visible within the result window after performing the test, the result is considered invalid. Instructions may not have been followed correctly. It is recommended to read the IFU again before re-testing the specimen with a new test device.

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1. Rothan HA, Byrareddy SN. The epidemiology and pathogenesis of coronavirus disease (COVID-19) outbreak. *J Autoimmun.* 2020; Feb 26:102433. doi: 10.1016/j.jaut.2020.102433.
2. Guo YR, Cao QD, Hong ZS, et al. The origin, transmission and clinical therapies on coronavirus disease 2019 (COVID-19) outbreak-an update on the status. *Mil Med Res.* 2020; Mar 13; 7(1):11. doi:10.1186/s40779-020-00240-0.
3. Lai CC, Shih TP, Ko WC, et al. Severe acute respiratory syndrome coronavirus 2 (SARS-Cov-2) and coronavirus disease-2019 (COVID-19): The epidemic and the challenges. *Int J Antimicrob Agents.* 2020; Mar 55(3): 105924. doi: 10.1016/j.ijantimicag.2020.105924.
4. *In Vitro Diagnostic Assays for COVID-19: Recent Advances and Emerging Trends* (Sandeep Kumar Vashist, 2020 April 05: diagnostics)
5. Nano Research for COVID-19 (<https://dx.doi.org/10.1021/acsnano.0c02540>)
6. Coronavirus (COVID-19) Update: FDA Authorizes First Antigen Test to Help in the Rapid Detection of the Virus that Causes COVID-19 in Patients (Stephen M. Hahn M.D. 2020 May 09: Commissioner of Food and Drugs)
7. Current Intelligence Bulletin 13: Explosive Ande Hazard DHHS (NIOSH) Publication Number 78-127 August 16, 1976
8. CDC. Discontinuation of Transmission-Based Precautions and Disposition of Patients with COVID-19 in Healthcare Settings (Interim Guidance). (2020).
9. CDC. Duration of Isolation and Precautions for Adults with COVID-19. (2020).
10. Bullard, et al. Predicting Infectious Severe Acute Respiratory Syndrome Coronavirus 2 From Diagnostic Samples. *CID.* 2020; Nov 15;71 (10); DOI:10.1093/cid/ciaa638

**GLOSSARY OF SYMBOLS /
SYMBOLVERZEICHNIS / GLOSARIO DE
SÍMBOLOS / GLOSSAIRE DES SYMBOLES /
GLOSSARIO DEI SIMBOLI / GLOSSÁRIO DE
SÍMBOLOS / ГЛОССАРИЙ ОБОЗНАЧЕНИЙ**

	Temperature limitation Temperaturbegrenzung Limitación de temperatura Limitation de température Limitazione di temperatura Limitação de temperatura Температурный диапазон
	For in vitro diagnostic use only Medizinprodukt für in-vitro Diagnostik Sólo para uso diagnóstico in vitro Pour un usage de diagnostic in vitro uniquement Ad uso esclusivo diagnostico in vitro Somente para uso para diagnóstico in vitro Медицинское изделие для диагностики In Vitro
	Do not reuse Nicht wiederverwenden No reutilizar Ne pas réutiliser Non riutilizzare Não reutilizar Запрет на повторное использование
	Do not use if package is damaged Bei beschädigter Verpackung nicht verwenden No lo use si el paquete está dañado Ne pas utiliser si le colis est endommagé Non utilizzare se la confezione è danneggiata Não use se o pacote estiver danificado Не использовать при повреждении упаковки

	Lot Number Chargencode Número de lote Numéro de lot Numero di lotto Número de Lote Номер серии
	Catalog Number Artikelnummer Número de catalogo Numéro de catalogue Numero di catalogo Número no Catálogo Номер изделия по каталогу
	Consult instructions for use Gebrauchsanleitung beachten Consultar instrucciones de uso Consulter les instructions d'utilisation Consultare le istruzioni per l'uso Consultar as instruções de uso Обратитесь к инструкции по применению
	Keep dry Trocken aufbewahren Mantener seco Garder au sec Mantenere asciutto Manter seco Беречь от влаги
	Biological Risks Biologisches Risiko Riesgos biológicos Risques biologiques Rischi biologici Riscos Biológicos Биологическая опасность
	Use By Verwendbar bis Usar por Utiliser par Utilizzare per Usar até Использовать до

	Manufacturer Hersteller Fabricante Fabricant Produttore Fabricante Производитель
	Date of manufacture Herstellungsdatum Fecha de manufactura Date de fabrication Data di produzione Data de fabrica��o Дата производства
	Keep away from sunlight Von Sonnenlicht fernhalten Mantener alejado de la luz solar Tenir � l'�cart de la lumi�re du soleil Tenere lontano dalla luce solare Manter longe da luz solar Не допускать воздействия солнечного света
	CE mark CE Zeichen Marca CE Marquage CE Marcatura CE Marca CE Знак соответствия продукции техническим регламентам ЕС
	Contains sufficient for X tests Ausreichend f�r X Pr�fungen Contiene suficiente para X pruebas Contient suffisamment pour les tests X Contenuto sufficiente per X test Cont�m suficiente para X testes Содержимого достаточно для проведения n-количества тестов
	Caution Achtung Precauci�n Attention Attenzione Cuidado Осторожно! Обратитесь к инструкции по применению

STERILE EO	Sterilized using ethylene oxide Sterilisiert mit Ethylenoxid Esterilizado con óxido de etileno. Stérilisé à l'oxyde d'éthylène Sterilizzato con ossido di etilene Esterilizado com óxido de etileno Стерилизовано оксидом этилена
STERILE R	Sterilized using irradiation Sterilisiert durch Bestrahlung Esterilizado mediante irradiación. Stérilisé par irradiation Sterilizzato con irradiazione Esterilizado por irradiação Радиационная стерилизация
	Do not re-sterilize Nicht sterilisieren No volver a esterilizar Ne pas re-stériliser Non ristertilizzare Não reesterilize Не стерилизовать повторно
CONTROL -	Negative control Negative Kontrolle Control negativo Contrôle négatif Controllo negativo Controle negativo Отрицательный контроль
CONTROL +	Positive control Positive Kontrolle Control positivo Contrôle positif Controllo positivo Controle positivo Положительный контроль

TECHNICAL SUPPORT / TECHNISCHER SUPPORT / SOPORTE TÉCNICO / SUPPORT TECHNIQUE / SUPPORTO TECNICO / SUPORTE TÉCNICO / ОТДЕЛ ТЕХНИЧЕСКОЙ ПОДДЕРЖКИ	
Europe & Middle East / Europa & Mittlerer Osten / Europa y Medio Oriente / Europe & Moyen Orient / Europe & Middle East / Europa e Oriente Médio / Европа, Ближний Восток	+44 161 483 9032 EME.techsupport@abbott.com
Africa / Afrika / África / Afrique / Africa / África / Африка	+27 10 500 9700 ARCIS.techsupport@abbott.com
Russia & CIS / Russland & CIS / Rusia y CEI / Russie & CIS / Russia & CIS / Rússia e CEI / Россия и СНГ	+44 161 483 9032 ARCIS.techsupport@abbott.com
Asia Pacific / Asien Pazifik / Asia Pacífico / Asie Pacifique / Asia Pacific / Ásia-Pacífico / Азия, Океания	+61 7 3363 7100 AP.techsupport@abbott.com
Latin America / Lateinamerika / America Latina / Amérique Latine / Latin America / América Latina / Латинская Америка	+57 1 482 4033 LA.techsupport@abbott.com
Canada/Kanada/Canada/ Canada/Canadá/Canadá/ Канада	+1 800 818 8335 Canproductsupport@abbott.com

 **Abbott Rapid Diagnostics Jena GmbH**
Ortenweg 1, D-07743 Jena, Germany
abbott.com/poet

Date Issued : 2021.03
41FK11/41FK21-07-A3

© 2021 Abbott. All rights reserved.

All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

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