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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality piece of equipment from KARL STORZ.

This manual is intended to serve as an aid in the proper installation, connection, and operation of the TELECAM C3. All the necessary details and actions are clearly explained. We therefore ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for easy reference.

WARNING: KARL STORZ medical

devices are not cleaned and disinfected when delivered. These hygiene measures must be performed before the first use/set-up. Cleaning/disinfection of the surfaces is part of the ongoing disinfection measures.

KARL STORZ is continuously working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, figures and descriptions in this manual.

- © All product illustrations, product descriptions and texts are the intellectual property of KARL STORZ GmbH & Co. KG.
Their use and reproduction by third parties require the express approval of KARL STORZ SE & Co. KG.
All rights reserved.

Rx only

Federal (USA) law restricts this device to sale by or on the order of a physician.

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It is recommended that the suitability of the products for the intended procedure be checked prior to use.

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Controls, displays, connectors and their uses

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2 Controls, displays, connectors and their uses

①	On/Off button
②	White balance
③	Camera connection C-Line
④	Camera connection X-Line
⑤	USB 2.0 (1x)
⑥	Potential equalization
⑦	DVI Out (2x)
⑧	USB 2.0 (4x)
⑨	LAN interface
⑩	Power cord socket

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3 Symbols employed

	Follow instructions for use
	Standby/On/Off
	Type CF/Defib. resistant applied part according to IEC 60601-1
	Type BF applied part according to IEC 60601-1
	White balance
	USB 2.0
	Potential equalization connector
	Control of pollution caused by electronic information products (China RoHS)
	This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	In accordance with US federal law (21 CFR 801.109), this product may only be sold to or on prescription from a licensed physician.
	Caution
	Manufacturer
	Catalogue number
	Serial number
	Number of products in the product packaging

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7.2.3	Function space	22
7.2.4	Setup menu	22
7.2.5	Information area	24
7.3	Quick menu configuration	26
7.4	Access and security concept	28
7.4.1	Encryption of sensitive menus via password	28
7.4.2	Role-based access system via individual user accounts	28
7.4.3	Memory functionality	30
7.5	Patient data security	31
7.5.1	Patient handling	32
7.5.2	Labeling	33
7.5.3	Export mode	33
8	Service and repair	34
8.1	Reprocessing	34
8.1.1	General safety instructions	34
8.1.2	Wipe-down disinfection of devices	35
8.2	Maintenance and safety check	36
8.2.1	Maintenance	36
8.2.2	Safety check	36
8.3	Servicing and repair	37
8.4	Disposal	37
8.5	Repairing devices	38
8.6	Infection prevention	38
8.7	Limitation of liability	39
8.8	Warranty	39
9	Technical description	40
9.1	Technical data	40
9.2	Standard compliance	41
9.3	Directive compliance	41
9.4	Software ownership and licensing	41
10	Spare parts, recommended accessories	43
10.1	Compatible camera heads for X-Line connection	43
10.2	Compatible camera heads for C-Line connection	43

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4 Description of the device

The TELECAM C3 is a camera control unit for endoscopic image display and documentation. A new, intuitive user interface assures simple and quick access as well as rapid understanding of menu navigation.

Functions are described in detail in chapter 7.

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5 Safety instructions

5.1 Explanation of Warnings and Cautions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

Please read this manual and follow the instructions carefully. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system. To make these words stand out more clearly, they are accompanied by a pictogram.

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 **WARNING:** Consult and follow the storage, transport and operating conditions prior to storing, transporting or operating, respectively. If the device is stored in a warm, humid environment, ensure that it has suitably equilibrated to the operating environment prior to use.

 **WARNING:** To reduce the risk of electrical shock, do not remove cover of unit. Refer servicing to qualified personnel. Removal of covers by unauthorized personnel will void the equipment's warranty.

 **WARNING:** Perceived electrical shocks may be experienced when using the camera system in conjunction with electrosurgical units.

 **WARNING:** High energy radiated light through endoscopes may give rise to high temperatures in front of the light outlet and to the tip of the endoscope. To minimize the risk of burns, always adjust the light source to the minimum illumination intensity necessary to achieve optimum illumination of the endoscopic scene when coupled to the video camera.

 **WARNING:** Refer to the appropriate section of this manual for validated cleaning and disinfection instructions.

 **CAUTION:** Do not store liquids on or above the unit.

 **NOTE:** Disposal of these system products, at the end of their useful life, must be in accordance with local regulations.

 **NOTE:** Do not discard as unsorted municipal waste.

 **NOTE:** Discard as electrical/electronic waste; recycle or reuse accordingly.

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 **WARNING:** Under some Single Fault Conditions per the IEC 60601-1 standard, the temperature of the equipment surfaces could get hot and there is a possible risk of a burn if touched.

 **WARNING:** No modification of this equipment is allowed.

 **NOTE:** Consult local authorities for reuse/recycle instructions.

 **NOTE:** The presence of maximum mains voltage on SIP/SOPs is not considered a normal condition.

 **NOTE:** Additional equipment connected to medical electrical equipment must be able to prove their compliance with the respective IEC or ISO standards (e.g., IEC 60950 for data processing equipment). Furthermore, all configurations shall comply with the standard requirements for medical systems (see Clause 16 of the 3rd Ed. of IEC 60601-1). Anybody connecting additional equipment to medical electrical equipment is a system configurator and is therefore responsible for the system's compliance with the standard requirements for systems. Please note that local laws take priority over the above-mentioned standard requirements. If in doubt, please consult your local specialist dealer or the technical service department. (Standard / directive references: IEC 60601-1+A1+A2: 6.8.2.c, 19.2.b, 19.2.c, IEC 60601-1: 7.9.2.5, 8.1, 16.2.d, MDD 93/42/EEC: Annex I clause 13.6.c).

 **NOTE:** For use with the product, KARL STORZ recommends the accessories listed in chapter 11. If third-party accessories are used, the user must check the correct operation himself.

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-  **WARNING:** To ensure safe operation, do not simultaneously touch the device output connectors and the patient.
-  **WARNING:** When endoscopes are used with energized endoscopically used accessories, the patient leakage currents may be additive. This is particularly relevant when using type CF devices. In this case, a Type CF endoscopically used accessory should be used in order to minimize the total patient leakage current.
-  **WARNING:** During discharge of a defibrillator, the image displayed may take up to 4 seconds to recover and in some cases a power cycle may be required to restore the image.
-  **WARNING:** The device should be placed on a flat surface or in a video cart. Do not block air inlets/outlets as overheating may result in fire hazard or device failure during procedure.
-  **WARNING:** Failure to perform proper white balance of system may result in color distortion.
-  **WARNING:** The operator is responsible to limit physical access to the TELECAM C3 from unauthorized persons.
-  **WARNING:** In the case of loss of surgical image for greater than 3 seconds, a power cycle may be required to restore the image.
-  **WARNING:** When using high frequency surgical equipment, keep the working part of the active electrode in the field of view of the operator to avoid accidental high frequency burns. Avoid contact between the active electrode and metal parts of the endoscope.

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CAUTION: Regular back-ups must be performed to safeguard data. KARL STORZ is not liable for loss of data. Responsibility for the data lies with the operator.

NOTE: All data stored internally on the device is stored in encrypted form. The data must be exported onto a USB storage medium before servicing.

NOTE: Upon saving the data to USB, patient data which is legible as plain text is stored on the storage medium.

5.3 Normal use

5.3.1 Intended use

The TELECAM C3 is a camera control unit for use with camera heads or videoscopes and is used for visualizing and documenting endoscopic and microscopic procedures.

The TELECAM C3 has no direct contact with the human body.

5.3.2 Contraindications

No contraindications relating directly to the medical device are currently known. The responsible physician must decide whether the foreseen application is admissible based on the general condition of the patient.

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5.3.3 Patient profile

The use of the TELECAM C3 is not limited to a certain patient profile (sex, age, weight etc.). The attending physician must assess a patient's health and declare it adequate.

5.3.4 User profile (specialist physician)

- Recognized medical skills for endoscopic interventions in the relevant specialization
- Adequate powers of comprehension to assess the situation rationally at all times during the application
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- No physical impairments that could diminish perception of activation and visual message

5.3.5 User profile (assistant)

- Completed training as medical assistant
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- No physical impairments that could diminish perception of activation and visual message

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5.4 Safety precautions at the site of installation

The device may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations. It is not intended for use in potentially explosive atmospheres. This means, among other things:

- Wherever easily combustible and explosive inhalation anesthetics or mixtures thereof are used, the device must not be operated in the demarcated hazard zone.
- This also applies for easily combustible and explosive chemicals, e.g., skin disinfectants and fast-acting surface disinfectants.

The device is equipped with a connector for potential equalization. It should be connected in accordance with the applicable national regulations.

WARNING: The device should be placed on a flat surface, in a video cart or fixed in place. Do not block air inlets/outlets as overheating may result in fire hazard or device failure during procedure.

WARNING: Ensure sufficient air circulation for the device.

WARNING: The device is not intended for use in potentially explosive atmospheres. Do not operate the device within the demarcated hazard zone while explosive anesthetic gases are in use.

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5.5 Safety precautions when operating the device

It is the user's responsibility to make sure the device is safe and operating properly before use.

NOTE: The documentation accompanying the accessories used must also be observed.

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6 Initial commissioning

6.1 Unpacking

Carefully remove the TELECAM C3 and accessories from their packaging. Check the delivery for missing items and evidence of shipping damage.

Please file any complaints with the manufacturer or supplier immediately. If possible, retain the original packaging materials for later use; these can come in handy if the device has to be transported.

6.2 Basic equipment

- 1 TELECAM C3 TC100
- 1 Power cord
- 1 instruction manual

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1. Place the device on a level surface. Do not cover the lateral air vents.

WARNING: The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

6.3.2 Connecting the potential equalization

NOTE: The device's ground line should be installed by a qualified electrician.

The TELECAM C3 is equipped with a connector for potential equalization.

1. Connect the potential equalization line with the connector for potential equalization on the device.
2. Connect the potential equalization line to the connector of the treatment room.

6.3.3 Connecting the monitor

1. Connect the monitor to the TELECAM C3 via one of the DVI outputs on the rear of the device. The primary working screen should be connected to DVI1.

NOTE: The connected monitor must support a resolution of 1920 x 1080 pixels.

6.3.4 Connecting USB devices

The USB interfaces of the TELECAM C3 are suitable for connecting up the following peripherals:

- Keyboard/mouse
- Printer
- Footswitch
- USB stick, FAT32 formatted
- USB to ACC adaptor

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6.3.5 Connecting the power supply

! WARNING: Only operate the device with the voltage stated on the manufacturer's identification plate.

1. Before connecting the devices to the electrical supply line, check that the values on the power sockets concur with the values on the device's identification plate.
2. Insert the device plug of the power cord into the power supply cord receptacle of the TELECAM C3.
3. Insert the TELECAM C3 power plug into the power sockets.

! WARNING: Always use a hospital grade power cord with this camera.

6.3.6 Switching on

To switch on the device, press the ON/OFF key on the front of the device.

6.3.7 Connecting the camera head/ videoscopes

Insert the camera connector into the corresponding socket.

i NOTE: Please observe the compatibility list (chapter 10)

i NOTE: Please note that 2 cameras can be connected to the device at the same time; however, only one camera can be displayed. If you have connected 2 cameras, you can select the relevant image with the 'Swap Cameras' button in the menu (see chapter 7.2.2).

i NOTE: Please ensure that the camera connector is always attached properly and securely before commencing the intervention.

! WARNING: Always ensure the connector is completely clean and free of debris. If moisture is present, dry thoroughly with sterile towel prior to insertion.

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6.3.8 White balance

NOTE: Perform a white balance before each use in order to adjust the perception of white by the camera and the ability to display colors properly.

1. Ensure that the imaging instrument* is connected to the TELECAM C3 and displays a live image on the monitor.
2. Connect the light source to the imaging instrument via the light cable.
3. Point the tip of the imaging instrument at any clean white surface (e.g., white gauze).
4. Make sure that the light source emits enough light to illuminate the white surface sufficiently. The monitor should display a completely white and well illuminated live image.
5. Press the  button on the front of the TELECAM C3 or the  button in the live menu to perform the white balance.
6. Look at the monitor to determine whether the white balance was successful or not.

*A camera head with a connected endoscope, videoscope or exoscope

6.3.9 Function test

1. Check the TELECAM C3 and the camera head or videoscope for external damage.
2. Ensure that the camera head cable and/or the videoscope cable is not damaged or kinked.
3. Point the camera head and/or the videoscope towards an object and check the quality of the display on the monitor.

NOTE: The device and accessories may only be used once the test for proper functioning has been successfully completed.

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7 Operating instructions

7.1 Basic information on operation

Operation of the TELECAM C3 is divided into the following 4 areas:

- Patient area
- Quick menu
- Function space and
- Setup

In addition, you will find an information area (information bar or information panel) in the top right-hand corner.

7.1.1 Camera head buttons

All compatible camera heads and all X-Line videoscopes can be used to call up the programmed functions and for menu control.

- ② Scroll up: Scrolls up the menus, options.
- ③ Menu/selection: Access to the quick menu, with this button (as when pressing the ENTER button) you also activate a marked option on the computer keyboard.
- ④ Scroll down: Scrolls down the menus, options.

7.1.2 Keyboard and mouse

The keyboard and mouse can be used to navigate the menus. A keyboard can be used to enter data.

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7.2.2 Quick menu

The quick menu allows you to choose between different camera functions and change the settings of the device. Functions which are not available are shown in gray.

If you operate the device using the camera head buttons or the keyboard, you access the quick menu by pressing the enter button on the keyboard or by pressing the menu button on the camera head/videoscope.

NOTE: You can change the configuration of the icons in the quick menu via the setup menu (see chapter 7.3).

NOTE: Some icons in the quick menu are shown in gray until a camera head is connected.

Description of functions of the symbols used

White balance
A white balance is performed here.

NOTE: Perform a white balance before every use.

Freeze
You can freeze the image here; during this time the live image is shown in the top right-hand corner of the monitor.

Zoom
This enables digital magnification of the display.

Alignment
The displayed image can be flipped vertically or horizontally and rotated 180° here.

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Enhance

The digital fiberscope filter (anti-moiré/anti-grid filter) can be set to two levels (A, B).

With these filters, interfering artifacts, such as the moiré effect and display of the fiberscope grid, are reduced. If the honeycombed grid structures are relatively small in the image, Filter A is recommended.

Filter B is particularly recommended in the case of larger honeycomb structures.

Light

The light source can be switched on and off here.

NOTE: If a videoscope with an integrated light source is used, the integrated LED of the videoscope can be switched on and off via this function.

Exit

You can exit the quick menu here. Alternatively, you can also exit the menu by clicking the 'Esc' button on the keyboard.

NOTE: Upon pressing the 'Exit' button, the system returns to its previous state.

Image capture

Still images can be captured here.

NOTE: To save images, a USB storage medium is recommended. At the same time, the images can be printed using a connected printer.

NOTE: In addition to or in place of a USB storage medium, the AIDA® documentation unit can be used to capture still images. These are saved on the AIDA® following the image capture.

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Video recording
Video recordings can be started and stopped here.

Training mode
A circle can be displayed in the centre of the image here. The training mode can be used specifically for endoscopy training with the TELECAM C3.

Camera brightness
The brightness of the camera can be set here.

Print
Print now: Immediately prints all images in the print queue.

NOTE: This button only appears if a printer is connected.

Swap cameras
Here you can switch between two connected cameras without having to disconnect one camera.

NOTE: This button only appears if two cameras are connected at the same time.

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7.2.3 Function space

Just like the quick menu, the function space (see chapter 7.2) allows you to choose between different camera functions and change the settings of the device (via the setup menu). It automatically contains all functions which were not configured in the quick menu.

NOTE: Some icons in the function space are shown in gray until a camera head is connected.

The function space is called up using the button shown.

In the setup menu (see 7.2.4) you can access the function space by pressing the setup button (below).

7.2.4 Setup menu

You can change the general TELECAM C3 settings in the setup menu.

Display language
The language of the user interface can be selected here.

Keyboard language
The keyboard language of the connected external keyboard can be selected here.

DVI signal output
You can select the image refresh rate between 50 Hz and 60 Hz for the DVI output here.

NOTE: After changing the frequency, the device must be restarted.

Date & time
The date, time and display format can be set here.

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Patient management
The settings for handling patient data can be managed here. Please see chapter 7.5 for details.

Module information
All the necessary information on the system and the connected camera heads and videoscopes can be called up here.

Quick menu settings
The quick menu can be configured here (see chapter 7.3).

SW update
The device software and the connected camera heads and videoscopes can be updated here.

Software licenses
The device contains open source software. A list of the open source software included in the product, and the open source software license terms can be viewed here.

Logged events
Audit logs
You can view the audit log here or export it via USB.

Factory settings
The device can be reset to these factory settings here, note the information in chapter 7.4.

Access & protection
You can manage access to sensitive areas of the system here.
You will find more details in chapter 7.4 Access and security concept

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Interner Speicher

7.2.5 Information area

The information area in the top right-hand corner of the screen comprises the information bar and the information panel (see chapter 7.2).

The information bar offers a continuous overview of the system status. You will find more information on the system status and the connected accessories in the information panel. The information panel is folded down by pressing the 'i' symbol in the system bar (see figure).

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Overview of connected camera head/videoscope (incl. button assignment)

If a compatible camera head is connected to the TELECAM C3, you can call up the article number and name as well as the camera head button assignment here.

The camera head button assignment is shown true to original; the assignment can be altered using the 'Settings' button.

NOTE: The camera head button assignment cannot be changed on all videoscopes, note whether a fixed assignment is already printed on the videoscope.

Overview of a connected footswitch

If a compatible footswitch is connected to the TELECAM C3, you can call up the article number and name as well as the button assignment here.

Overview of storage medium used

If a storage medium is selected and connected to the TELECAM C3, you have the option here of calling up the currently used and free storage space.

Should you use the internal storage, you can transfer the recorded data to a USB storage medium here.

NOTE: Internally saved data is deleted after successful transfer to the USB storage medium in order to create space for new recordings.

You will find more information on the device storage options in chapter 7.4.3 Storage functionality

Overview of a connected printer

If a compatible USB printer is connected to the TELECAM C3, you have the option here of calling up the printer configuration and changing this if necessary. The images per page and the number of copies can be changed and it can also be set whether images are to be printed out automatically during an ongoing treatment.

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The configuration menu is divided into 3 columns. Here functions can be moved, replaced, added or removed from the quick menu.

NOTE: You can add up to 8 functions to the quick menu.

NOTE: The configuration is saved by restarting the device.

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- ① **Quick menu:**
The current functions are listed here (see 7.2.2).
-  Press this button to add a new function to the quick menu.
- ② **Command:**
The functions in the quick menu can be configured using the 'Command' buttons. To this end, select a function in the (left) column ①.
-  Press this button to move the selected function upwards in the quick menu.
-  Press the 'Delete' button to remove the selected function from the quick menu.
-  Press the 'Replace' button to select another function (see ③) in place of the current one.
-  Press this button to move the selected function downwards in the quick menu.
- ③ **Select Function:** The functions which are not currently configured in the quick menu are listed here.

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7.4 Access and security concept

The TELECAM C3 offers two options for protecting the device and sensitive menu areas from access from unauthorized parties. One of the two options must be selected when configuring the system in order to work with sensitive information such as patient data.

i NOTE: Please note that for security reasons, the device must be reset to the factory settings if the password is lost. As such, in addition to all settings, the patient data is also deleted from the device and the device is reset to the default patient privacy settings.

7.4.1 Encryption of sensitive menus via password

By selecting this option, certain menus are protected by a password. The standard password is [90290245] and must be changed.

7.4.2 Role-based access system via individual user accounts

Upon selecting this option, sensitive areas of the menu as well as device functions can only be selected by certain users. Authorization to call up these menus and functions is linked to the specific roles of the user accounts.

Configuration of the role-based access system

When you configure the role-based access system, you must firstly create an administrator account. Moreover, in the first step you must specify password rules which apply to all further passwords. Then, with the help of the standard password, you can create a new password and subsequently make the further settings.

i NOTE: The 'Administrator' role has full access to all device functions and settings.

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NOTE: In addition to the Administrator, the roles of 'OR Personal' and 'Technician' can also be selected. The 'OR Personal' role enables access to patient data, but not to the patient data security and memory areas. The role of 'Technician' enables access to device settings in order, for example, to perform updates without providing access to patient data or other sensitive areas.

NOTE: All safety-relevant actions are logged and recorded by the system. KARL STORZ assumes no legal liability or responsibility for the loss, accuracy, damage or disclosure of patient data which was managed using TELECAM C3.

1. Configuration of further user accounts
You can create further user accounts, edit existing accounts and also delete accounts under this menu item.
A name, corresponding role and a password must be specified for each user account.

NOTE: Limitation of the number of user accounts.

2. Login message
You can enter a customized message here which is shown in the login area before granting access to sensitive areas.

3. Password rules
Here you can set whether and which rules are to be applied when assigning passwords.

4. Automatic logout
You can set here whether the system should log out a user after expiry of an individually set period of inactivity.

5. Deactivation
Upon deactivating the role-based access system, the encryption of sensitive menus via password is automatically activated, i.e., an individual password must be created.

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Use of the role-based access system

After you have activated the role-based access system, a further menu item with the name 'Login' appears in the main screen next to the information area.

Here all system users can log in using their user account and password. If no users are logged in, certain areas such as the patient area and the settings for patient data security and memory cannot be called up.

7.4.3 Memory functionality

The TELECAM C3 comprises 2 storage options:

- Internal storage and
- USB storage

Depending on the configuration of the recording and storage options in the setup menu under the item Access and security, images and videos with patient data can be captured in the folder and file name. The type of storage can also be selected here.

For storage on USB, the internal storage is used for buffering. This means data is firstly buffered in the internal storage and then transferred to the USB.

NOTE: The transfer status of the data can be checked in the information area, see chapter 7.2.5 Symbols in the information area.

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In order to store recordings without any loss of data, the 'Complete treatment' button in the patient area must be pressed. The pending data is then transferred and saved properly.

If no USB storage medium is available, the internal storage of the TELECAM C3 can also be used. Approx. 50 GB of storage space are available for internal storage.

NOTE: Internally stored recordings can only be removed from the internal storage via USB export. Please note that if large volumes of data are saved in the internal storage, a correspondingly large USB storage medium must be available.

NOTE: The transfer of large data volumes to a USB storage medium can take some time. The current status of data transfer can be checked in the information area.

7.5 Patient data security

NOTE: The TELECAM C3 system does not contain patient data as standard. The patient data administration functions (e.g., entry, access, display, storage and transfer) must be activated by the user.

NOTE: Upon activating the patient data administration functions, the user assumes responsibility for all patient-related data which are managed using TELECAM C3 as per the applicable national data protection laws.

NOTE: KARL STORZ assumes no legal liability or responsibility for the loss, accuracy, damage or disclosure of patient data which is managed using TELECAM C3.

In the Patient Data Privacy settings you can configure all settings related to the patient information.

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Patient information of any kind is only available if the corresponding settings have been made.

The configuration of the 'Patient Data Privacy' menu is password-protected as standard. You will find further information on password protection in chapter 7.4. Access and security concept.

The 'Patient Data Privacy' menu is divided into 3 areas:

- Patient handling
- Labeling and
- Export mode

7.5.1 Patient handling

1) Worklist setup

- No Worklist / Patient Entry (selected as standard)
- Local Worklist

Here you can establish whether patient data may be used (this is not permitted as standard).

i NOTE: When you are working without entering patient data, the data in the info panel must be exported via 'Transfer to USB'.

2) Capture tools

- No Capture
- Capture Allowed (selected as standard)

Here you can specify whether and how images and videos are to be captured. The following options are available: no capture, capture of images and videos labeled with patient information and the anonymous capture of media files.

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3) Storage limit

- Infinite
- Power cycle
- 1 day
- 7 days (selected as standard)
- 30 days
- 90 days

Here you can specify how long data sets are to be stored on the device.

i NOTE: Please note that both patient data sets in the worklist and media files (recorded images and videos) are deleted if this is set here. The data can no longer be restored.

7.5.2 Labeling

Allow patient data on prints

- No patient data on prints (selected as standard)
- Patient data on prints

Here you can determine whether prints will be labeled with patient information or not.

7.5.3 Export mode

Selecting the Export mode

- No export
- USB (or USB anonymous) export (selected as standard)

Here you can establish whether data export to USB is to be permitted.

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8.1 Reprocessing

8.1.1 General safety instructions

 **WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

 **WARNING:** When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

 **WARNING:** When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

 **WARNING:** National laws and regulations must be observed.

 **CAUTION:** Do not allow water or any liquid to drip inside the unit(s).

 **CAUTION:** Damage to product. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be requested from hygiene@karlstorz.com.
*Not available for U.S. Customers

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**8.1.2 Wipe-down disinfection of
devices**

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry low lint cloth.

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8.2 Maintenance and safety check

8.2.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the device's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

8.2.2 Safety check

WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.

Visual inspection

- Check the device and accessories for any mechanical damage which may impair functionality.
- Check that the inscriptions relevant to safety are legible.

Electric measurements

- Protective grounding measured according to IEC 62353 (current edition), and checked for compliance with the limits specified in the standard.
- Touch Current, Earth Leakage Current, and Patient Leakage currents measured according to IEC 62353 (current edition), and checked for compliance with the limits specified in the standard.

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78579 Neuhausen

Servicehotline: +49 7461 708-980

E-Mail: technicalsupport@karlstorz.com

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8.5 Repairing devices

Only personnel authorized by the manufacturer may repair devices. Interventions described in these instructions for use are exempt from this rule.

In Germany, you can contact the following address for repairs:

KARL STORZ SE & Co. KG

Abt. Reparaturservice

Take-Off Gewerbepark 83

78579 Neuhausen, Germany

Service hotline: +49 7461 708-980

E-mail: technicalsupport@karlstorz.com

In other countries, the respective KARL STORZ branches or specialist dealers are responsible.

8.6 Infection prevention

Contaminated devices may not be shipped. To prevent contact infections and airborne infections, devices must first be decontaminated. KARL STORZ reserves the right to send back contaminated devices.

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8.7 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications or repairs have been performed by a person or persons duly authorized by KARL STORZ,
- all electrical installations at the location in which the device is connected and operated meet the applicable laws and standards, and
- the device has been used in accordance with its operating instructions at all times.

8.8 Warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see chapter 'Subsidiaries'), even during the warranty period.

Opening of the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair or modification performed during the guarantee period shall void all warranty.

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9 Technical description
9.1 Technical data

TELECAM C3	TC100
Line voltage	100...240 VAC
Power frequency	50/60 Hz
Power consumption	70 VA
Dimensions	305 mm x 77 mm x 339 mm (L x W x D)
Weight	2.75 kg
Storage/transport temperature	-20°C to +60°C (-4°F to 140°F)
Humidity (non-condensing)	Operation: 20% ... 95% Storage/transport: 5% ... 95%
Air pressure	Storage/transport: 500 hPa – 1080 hPa
Operating temperature	0°C to 40°C (32°F to 104°F)
Operating rated altitude NN	3,000 m
Technical data Memory device	
Memory interface	USB
Image format	JPEG
(Audio)/Video format	MPEG-4
Internal Storage Size	50 GB
Interfaces	
Service interface	RJ45
Video interface	2x DVI-D output
Footswitch connection	USB

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- Type of protection against electric shock:
Protection Class I
- Protection of camera applied part against electric shock:

Depending on the camera connection used,
type BF or CF, defibrillation-proof

- Type of protection against moisture: IPX0

① NOTE: Please observe the electromagnetic compatibility information in chapter 11.

9.3 Directive compliance (for TC100)

According to Medical Device Directive (MDD):

This medical device belongs to Class I
This medical device bears a CE mark in accordance with MDD 93/42/EEC.

9.4 Software ownership and licensing

With respect to KARL STORZ products containing software components, the customer has a nonexclusive, limited, non-transferable license to use the programmed logic, computer programs and/or software supplied by KARL STORZ, in connection with, and incorporated into, KARL STORZ products (Software). The user right is only valid for the delivered form and only for use in accordance with the printed instruction manual from KARL STORZ for the products sold to the customer (and for no other product). The software, and all modifications, enhancements and upgrades thereto, will, at all times, remain the property of KARL STORZ. The customer may not duplicate, copy, reverse-engineer, de-compile, or disassemble the software or in any way modify the software. The customer has no right to create derivatives of the software, and the customer may

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not copy, decipher or re-create the source code of the software. Any and all such modifications or enhancements to the software by the customer, in contravention of this license, will immediately become the sole property of KARL STORZ.

The customer hereby acknowledges and agrees (i) that the purchase of KARL STORZ products does not constitute a sale of the software, (ii) that the software is the property of KARL STORZ, (iii) that the customer neither owns nor acquires any claim or intellectual property right in or to the Software, and (iv) KARL STORZ retains all right, title, and interest in and to the software, at all times, regardless of the form or media in or on which the original or other copies of the software may exist. In the event of a failure of the customer, or its agents, employees or representatives, to comply with any terms and conditions of the License herein granted, the license will, without any further action by KARL STORZ or any other party, immediately terminate.

Attribution notice

This device contains open source software. You can view a list of the open source software included in the product, and the open source software license terms, on the TELECAM C3 device itself. Select the tab Software licenses in the setup menu.

To obtain a copy of the source code being made publicly available by KARL STORZ Imaging, Inc. («KARL STORZ») related to open source software used in this KARL STORZ product, you should send your request, along with the device and serial number, in writing by email to opensource@karlstorz.com or by postal mail to:
KARL STORZ Imaging, Inc.
Attention: Open Source Software Inquiries
1 S. Los Cameros
Goleta, CA 93117
USA

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10 Spare parts, recommended accessories

NOTE: For use with the product, KARL STORZ recommends the accessories listed in this chapter. If third-party accessories are used, the user must check the correct operation himself.

NOTE: Not all accessories are available in all regions.

10.1 Compatible camera heads for X-Line connection

Article/model	Description
TH110	IMAGE1 S™ HX
TH111	IMAGE1 S™ HX-P
TH115	IMAGE1 S™ D1
TH116	IMAGE1 S™ D1

10.2 Compatible camera heads for C-Line connection

Article/model	Description
TH130	H1

10.3 Compatible videoscopes for X-Line connection

Article/model	Description
ENT	
11101VP/VN	CCD Video Rhino-Laryngoscope
11101VPS/VNS	Strobo Video Rhino-Laryngoscope
11101HD	HD Video Rhino-Laryngoscope

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**Spare parts,
recommended accessories**

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UROLOGY	
11278VSU/VSUE	Video Uretero-Renoscope FLEX-XC
11278VS/VSE	Video Uretero-Renoscope FLEX-XC
11272VP/VN	Flexible Video Urethro-Cystoscope
11272VPU/VNU	Flexible Video Urethro-Cystoscope
11272VH/VHU	Flexible Video Cystoscope HD-VIEW
PROCTOLOGY	
13912PKS/NKS*	TROIDL Flexible Rectoscope
SURGERY	
11292VP/VPU	Flexible Video Choledochoscope
11292VS/VSU	Flexible Video Choledochoscope (new sensor)
VETERINARY	
60714PKS/NKS*	Veterinary Video Endoscope
60914PKS/NKS*	Veterinary Video Endoscope
60118PKS/NKS*	Veterinary Video Endoscope
60125PKS/NKS*	Veterinary Video Endoscope
60130PKS/NKS*	Veterinary Video Endoscope
60332PKS/NKS*	Veterinary Video Endoscope
60511PKS/NKS*	Veterinary Video Endoscope
60278VS/VSU*	FLEX-XC Veterinary

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① **NOTE:** Not all accessories are available in all regions.

**10.4 Compatible videoscopes for
C-Line connection**

Article/model	Description
ANESTHESIA	
8403XS	C-MAC® S Imager
8403XSI	C-MAC® S Imager
8403XSP	C-MAC® S Pediatric Imager
8403EXC	C-MAC® Video Laryngoscope MAC #0
8403KXC	C-MAC® Video Laryngoscope MAC #2
8403AX	C-MAC® Video Laryngoscope MAC #3
8403AXC	C-MAC® Video Laryngoscope MAC #3
8403BX	C-MAC® Video Laryngoscope MAC #4
8403BXC	C-MAC® Video Laryngoscope MAC #4
8403HX	C-MAC® Video Laryngoscope D-BLADE
8403HXP	C-MAC® Video Laryngoscope D-BLADE PED.
8403DXC	C-MAC® Video Laryngoscope MIL #0
8403GXC	C-MAC® Video Laryngoscope MIL #1
8403MXC	C-MAC® Video Laryngoscope MIL #2
11301ABX	Flex. Intubation Videoscope
11302BDX	Flex. intubation videoscope (as of serial number 60002)
11303BNX	Flex. Intubation Videoscope
11304BCX	Flex. Intubation Videoscope
10331BX	C-MAC® VS

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i NOTE: Not all accessories are available in all regions.

ENT	
11102CM	CMOS Video Rhino-Laryngoscope
13303E	Video Esophagoscope
UROLOGY	
11272VE/VUE	CMOS video cysto-urethroscope

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Electromagnetic Compatibility (EMC) information

11 Electromagnetic Compatibility (EMC) information

The device is suitable for use in professional healthcare facilities.

Professional healthcare facilities include physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the high-energy shielded room of an ME system for magnetic resonance imaging).

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

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Accessories and lines for EMC compliance				
Type	Shielded	Length [m]	Ferrite	Use
Power cord	No	3	No	Power supply
Dual-pedal footswitch	Yes	2	No	Unit control, socket rear
USB keyboard with touchpad	Yes	2	No	Data entry
Camera	Yes	3	No	Image transmission
DVI cable	Yes	3	Yes	Signal transmission

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WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

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WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

NOTE: The emission characteristics of this device make it suitable for use in industrial areas as well as in hospitals (CISPR 11 Class A) and other professional healthcare environments. If the device is used in a domestic environment (for which CISPR 11 Class B is normally required), the device may not offer sufficient protection for radio transmission operation. The user might need to take mitigation measures, such as relocating or re-orienting the device.

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Table 1 – Compliance level for immunity tests
Guidelines and manufacturer's declaration – electromagnetic immunity

The device is designed for operation in the electromagnetic environment defined below. The user of the device should make sure that it is used in such an environment.

Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment – guidelines
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be made of wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity must be at least 30%.
Electrical fast transient / burst acc. to IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines 100 kHz repetition	± 2 kV/1 kV ± 1 kV for input and output lines 100 kHz repetition	The power supply quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV voltage outer conductor – outer conductor ± 2 kV voltage outer conductor – ground	± 1 kV voltage outer conductor – outer conductor ± 2 kV voltage outer conductor – ground	The power supply quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage dip:</u> Dip to 0% for 1 cycle @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 periods	<u>Voltage dip:</u> Dip to 0% for 1 cycle @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 periods	The power supply quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50 Hz/60 Hz	30 A/m at 50 Hz/60 Hz	If image distortion occurs, it may be necessary to position the device further from sources of electromagnetic fields or to install magnetic shielding. Before installing the device, the electromagnetic field should be measured to ensure that it is sufficiently low.
Interference to immunity test acc. to IEC 61000-4-3 for high-frequency electromagnetic fields	3 V/m 80 MHz to 2.7 GHz *Refer to table 2 for wireless Proximity RF field test levels	3 V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by high-frequency fields acc. to IEC 61000-4-6	3 Vrms to 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 Vrms in ISM band	3 Vrms to 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 Vrms in ISM band	

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Table 2

Test levels for Proximity fields from RF wireless communications equipment

Test frequency MHz	Frequency band MHz	Radio service	Modulation	Immunity test level V/m	Compliance level V/m
385	380 – 390	TETRA 400	Pulse modulation 18 Hz	27	27
450	430 – 470	GMRS 460 FRS 460	FM $\pm$ 5 kHz deviation 1 kHz sine wave	28	28
710	704 – 787	LTE band 13 & 17	Pulse modulation 217 Hz	9	9
745					
780					
810	800 – 960	GSM 800/900 TETRA 800 iDEN 820 CDMA 850 LTE band 5	Pulse modulation 18 Hz	28	28
870					
930					
1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	28	28
1845					
1970					
2450	2400 – 2570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE band 7	Pulse modulation 217 Hz	28	28
5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9	9
5500					
5785					

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Table 3 – Test levels for Radiated and Conducted Immunity Tests

Guidelines and manufacturer's declaration – electromagnetic immunity

The device is designed for operation in the electromagnetic environment defined below. The user of the device should make sure that it is used in such an environment.

Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment – guidelines
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{me} 150 kHz to 80 MHz	3 V _{me}	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ Where P is the rated power of the transmitter in watts [W] according to the information provided by the transmitter manufacturer and d is the recommended separation distance in meters [m]. Field strength of fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . $d = 1.2 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	Interferences may occur in the vicinity of equipment marked with the following symbol: 

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic waves is affected by absorptions and reflections of buildings, objects, and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the device is used exceeds the above compliance levels, the device should be monitored to ensure proper function. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 – Emission class and group

Guidelines and manufacturer's declaration – electromagnetic emissions

The device is designed for operation in an environment as defined below. The customer or user of the device should ensure that it is used in such an environment.

RF emissions acc. to CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions acc. to CISPR 11	Class A	The device is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Emission of harmonic oscillations acc. to IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions as per IEC 61000-3-3	complies	

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Table 5

Recommended separation distances between portable and mobile RF communications equipment and the device

The device is intended for use in an electromagnetic environment in which the RF disturbances are controlled. The customer or user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters whose maximum rated power is not listed in the table above, the recommended separation distance d in meters (m) can be estimated using the equation from the respective column, whereby P is the maximum rated power of the transmitter in watts (W) according to the transmitter manufacturer.

Note: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

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Compatibility with surgical HF instruments

This device was tested for compatibility with HF surgical devices. It has been tested according to IEC 60601-2-2 Annex BB.

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KARL STORZ Endoscopy Canada Ltd.
1717 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksea.com

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TOO KARL STORZ Endoscopy Kazakhstan
Saryarka, 6, BC "Arman", off. 910
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai – United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Interchange 21 Tower, Level 33
399 Sukhumvit Road,
North Klongtoey, Wattana,
10110 Bangkok, Thailand
Phone: +66 2 660 3669
E-Mail: info-th@karlstorz.com

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, MP Plaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd.
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Janguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch
Room 701A Building 5 & Room 501 Building 7,
No. 3000 Longdong Avenue, Pilot Free Trade Zone,
201203, Shanghai, People's Republic of China
Phone: +86 21 60339888, Fax: +86 21 60339808
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
Room 803-805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
610041, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhua District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch
Room 02B & 03 & 04A, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 87321281, Fax: +86 20 87321286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Indonesia
Sinarmas MSG Tower Level 37
Jl. Jend. Sudirman No. Kav. 21
Jakarta Selatan
DKJ Jakarta 12920
E-Mail: infoindonesia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info-jp@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
68 Waterloo Road, Macquarie Park NSW 2113
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (02) 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (02) 9420 0695
E-Mail: karlstorz@karlstorz.com.au

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Region "UROLOGY"

11272VE / VUE Videocystou

Compatible Disposable Video Endoscopes for C-Line Connector

09136101 FIVE S 3.5x65, sterile, single use
 09137001 Video esophagoscope CMOS SSU
 09133001 Video rhino-laryngoscope CMOS SSU

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
Figure 1. Example of nesting cardboard inserts

Table 2.

Characteristics	Value
Overall dimensions of a cardboard box, mm	428,6±43x 336,6±33x 161,9±16 (L x W x H)
Cardboard box thickness, mm	5

Table 3. Packing materials

Name of the component of the package	Brand and material
INSERT, BOX, SHIP	TC100 275 B-Flute Kraft (fine corrugated cardboard B)
SLEEVE, CHIP	TC100 SBS 24 (styrene-butadiene-styrene)
FILLER, PAD	TC100 200 B-Flute Kraft (fine corrugated cardboard B)
BOX, ACCESSORY	TC100 200 B-Flute Kraft (fine corrugated cardboard B)

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	B)
BOX, SHIPPER	TC100 200 C-Flute Kraft (medium corrugated cardboard with)
END CAP	TC100 200 B-Flute Kraft (fine corrugated cardboard B)

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Precautions

Explanation of Warnings and Cautions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system. Please read this manual and follow the instructions carefully. The words Warning, Caution and Note convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system. To make these words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician is at risk.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

NOTE: A Note indicates special information about operating the device or clarifies important information.

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Safety instructions

WARNING: Read this instruction manual thoroughly before using the equipment. Read the chapter on safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

WARNING: Risk of injury and damage to the products: Failure to observe and follow this instruction manual and all of the instruction manuals of the products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals carefully and always follow the instructions given precisely. Check the functioning of the products used in combination.

WARNING: Test this device prior to each surgical procedure to ensure that it functions correctly. In the event that the image becomes unusable during surgery, the camera may be disengaged from the endoscope and the procedure continued optically. If this is not possible, it is up to the discretion of the surgeon how best to proceed. Availability of a spare system is recommended.

WARNING: Danger of an electric shock! This device may only be connected to a socket with protective earthing. To reduce the risk of electric shock, keep the device out of reach of patients.

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WARNING: Grounding reliability can only be achieved when the equipment is connected to a 'Hospital Only' or 'Hospital Grade' outlet (i.e., approved for use in an operating room environment). The plug and cord should be routinely inspected. Do not use if damaged.

WARNING: Before each use or after a change of viewing modes/settings, the operator should check to ensure the view observed through the endoscope provides a live image (rather than a stored one) and has the correct image orientation

WARNING: Consult and follow the storage, transport and operating conditions prior to storing, transporting or operating, respectively. If the device is stored in a warm, humid environment, ensure that it has suitably equilibrated to the operating environment prior to use.

WARNING: To reduce the risk of electrical shock, do not remove cover of unit. Refer servicing to qualified personnel. Removal of covers by unauthorized personnel will void the equipment's warranty.

WARNING: Perceived electrical shocks may be experienced when using the camera system in conjunction with electrosurgical units.

WARNING: High energy radiated light through endoscopes may give rise to high temperatures in front of the light outlet and to the tip of the endoscope. To minimize the risk of burns, always adjust the light source to the minimum illumination intensity necessary to achieve optimum illumination of the endoscopic scene when coupled to the video camera.

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WARNING: Refer to the appropriate section of this manual for validated cleaning and disinfection instructions.

CAUTION: Do not store liquids on or above the unit.

NOTE: Disposal of these system products, at the end of their useful life, must be in accordance with local regulations.

NOTE: Do not discard as unsorted municipal waste.

NOTE: Discard as electrical/electronic waste; recycle or reuse accordingly.

WARNING: Under some Single Fault Conditions per the IEC 60601-1 standard, the temperature of the equipment surfaces could get hot and there is a possible risk of a burn if touched.**WARNING:** No modification of this equipment is allowed.

NOTE: Consult local authorities for reuse/recycle instructions.

NOTE: The presence of maximum mains voltage on SIP/SOPs is not considered a normal condition.

NOTE: Additional equipment connected to medical electrical equipment must be able to prove their compliance with the respective IEC or ISO standards (e.g., IEC 60950 for data processing equipment). Furthermore, all configurations shall comply with the standard requirements for medical systems (see Clause 16 of the 3rd Ed. of IEC 60601-1). Anybody connecting additional equipment to medical electrical equipment is a system configurator and is therefore responsible for the system's compliance with the standard requirements

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for systems. Please note that local laws take priority over the above-mentioned standard requirements. If in doubt, please consult your local specialist dealer or the technical service department. (Standard / directive references: IEC 60601-1+A1+A2: 6.8.2.c, 19.2.b, 19.2.c, IEC 60601-1: 7.9.2.5, 8.1, 16.2.d, MDD 93/42/EEC: Annex I clause 13.6.c).

NOTE: For use with the product, KARL STORZ recommends the accessories listed in chapter 11. If third-party accessories are used, the user must check the correct operation himself.

WARNING: To ensure safe operation, do not simultaneously touch the device output connectors and the patient.

WARNING: When endoscopes are used with energized endoscopically used accessories, the patient leakage currents may be additive. This is particularly relevant when using type CF devices. In this case, a Type CF endoscopically used accessory should be used in order to minimize the total patient leakage current.

WARNING: During discharge of a defibrillator, the image displayed may take up to 4 seconds to recover and in some cases a power cycle may be required to restore the image.

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WARNING: The device should be placed on a flat surface or in a video cart. Do not block air inlets/outlets as overheating may result in fire hazard or device failure during procedure.

WARNING: Failure to perform proper white balance of system may result in color distortion. **WARNING:** The operator is responsible to limit physical access to the TELECAM C3 from unauthorized persons.

WARNING: In the case of loss of surgical image for greater than 3 seconds, a power cycle may be required to restore the image.

WARNING: When using high frequency surgical equipment, keep the working part of the active electrode in the field of view of the operator to avoid accidental high frequency burns. Avoid contact between the active electrode and metal parts of the endoscope.

WARNING: When using two camera heads/videoscopes, always make sure that you have chosen the right source for the current procedure.

WARNING: This medical device provides and manages sensitive information. Therefore, access to the device must be restricted to persons who are authorized to handle patient data.

CAUTION: Regular back-ups must be performed to safeguard data. KARL STORZ is not liable for loss of data. Responsibility for the data

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