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Instruction for use

Nasal spray MARIMER

«APPROVED»

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Position
Responsible Pharmacist/Quality Director

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*I hereby certify accuracy, correctness and
reliability of this document
translated into Russian*

CONTENTS

NAME OF MEDICAL DEVICE	3
INFORMATION ABOUT MANUFACTURER	3
DEVELOPER / MANUFACTURER	3
MANUFACTURING SITE OF MEDICAL DEVICE	3
INTENDED USE OF MEDICAL DEVICE	3
INDICATIONS FOR USE	3
CONTRAINDICATIONS	3
SIDE EFFECTS	4
CONDITIONS OF USE	4
DESCRIPTION, COMPOSITION AND SUPPLY PACKAGE OF MEDICAL DEVICE	4
POSODOLOGY AND METHOD OF ADMINISTRATION	4
TECHNICAL DESCRIPTION OF MEDICAL DEVICE	5
PRECAUTIONS FOR USE	5
WARNING	5
STERILIZATION METHODS	5
CLEANSING	5
OPERATIONAL, STORAGE AND TRANSPORTATION CONDITIONS	5
SHELF-LIFE	6
ENVIRONMENTAL REQUIREMENTS	6
LABELLING	6
PROCEDURE FOR DISPOSAL OR DESTRUCTION OF MEDICAL DEVICE	6
LIST OF INTERNATIONAL STANDARDS	7
LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN	7
INFORMATION ABOUT DRUG PRODUCTS, MEDICINAL PRODUCTS AND PHARMACEUTICAL SUBSTANCES	7
WARRANTY	7
RECLAMATION	7

NAME OF MEDICAL DEVICE

Nasal spray MARIMER, composed of:

- Nasal spray MARIMER, 100 mL. – 1 pce.
- Package leaflet – 1 pce.

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Tel./Fax +33 2 31 47 15 15/+33 2 31 43 95 00

MANUFACTURING SITE OF MEDICAL DEVICE

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INTENDED USE OF MEDICAL DEVICE

Nasal spray MARIMER is intended for penetration, purification, release and moisturizing of the nasal passages and sinuses, for the complex treatment and prophylaxis of ENT diseases, viral infections (ARVI, flu), allergic and vasomotor rhinitis.

INDICATIONS FOR USE

Nasal spray MARIMER is used for the children from 1 year old and adults:

- for daily cleaning of the nasal cavities, including for the maintenance of the protective properties of the mucous membrane under conditions of increased dryness or air pollution (air conditioning, central heating, etc.);
- to moisturize the dry mucous membrane of the nasal cavity;
- for the complex treatment and prophylaxis of ENT diseases, viral infections (ARVI, influenza), allergic and vasomotor rhinitis.

CONTRAINDICATIONS

Children under 1 year.

SIDE EFFECTS

Side effects when using the device in accordance with the package leaflet have not been identified.

CONDITIONS OF USE

Used in medical facilities, as well as at home. Refers to personal use devices. Carrying out procedures with this device at home does not require special training and special skills.

DESCRIPTION, COMPOSITION AND SUPPLY PACKAGE OF MEDICAL DEVICE

Nasal spray MARIMER is an isotonic solution of seawater in a pouch inside an aluminum bottle equipped with an actuator with a valve.

The solution is administered self-sustained by pressing the actuator. Microdiffusion spray system provides deep and soft irrigation of the nasal mucosa, which contributes to the effective action of trace elements. The solution can be sprayed in any position. The trace elements in the composition of isotonic seawater solution improve the function of the ciliary epithelium, which increases the resistance of the nasal mucosa to pathogens and viruses.

Nasal spray MARIMER is packed in a 100 ml aluminum bottle with an actuator and a cap, and then placed in an amount of 1 pce. in a carton along with a package leaflet.

Composition per 100 mL: seawater – 28 mL, purified water – up to 100 mL.

The bottle is filled with compressed air (or an inert propellant gas), which is not in contact with the solution.

pH 6,5-8,5

Osmolality 270 - 310 mOsm/kg

POSODOLOGY AND METHOD OF ADMINISTRATION

Device can be used daily for an unlimited period of time.

One spray into each nasal passage from one to four times a day on demand.

Nasal spray MARIMER can be used in any position: lying, sitting or standing (when using in a child, it is recommended to take into account the general condition of the child and its age characteristics).

- Tilt the head to one side.
- Insert actuator into the upper nasal passage and spray isotonic solution.
- Repeat the procedure with another nasal passage.
- Straighten the head, wipe the excess from the nose with a tissue or blow the nose.

After each use, rinse the actuator with warm soapy water, rinse and wipe.

TECHNICAL DESCRIPTION OF MEDICAL DEVICE

Nasal spray MARIMER is manufactured in compliance with international and national standards. Isotonic seawater solution is packaged in an individual package, which consists of a bag with a valve placed in an aluminum bottle, an actuator and a cap.
Bag with valve (type Power pouch) made of material: polyethylene.

By the biological effect on the patient, the materials used satisfy the requirements of GOST R 52770 and the standards of the GOST ISO 10993 series. The product is biologically safe.

Product weight without consumer packaging 145 g (± 15 g)

Droplet size while spraying	Dv(10)- not more than 30 μ m Dv(50) – not more than 65 μ m Dv(90 – not more than 135 μ m
Weight of single sprayed dose	2.621 $\pm$ 5%

The mass of a spray nozzle is 4.6g $\pm$ 0.023g

The dimension of the cap is

- Height: 60mm $\pm$ 3.3mm
- Diameter 44.4 $\pm$ 0.2mm
- Weight: 8.2g $\pm$ 0.041g

PRECAUTIONS FOR USE

- Attention! Bottle is under pressure!
- Store at the temperature from 10 to 50°C.
- Keep out of reach of children!

WARNING

Hermetically sealed container: may explode when heated. Keep away from heat, sparks, open flames, hot surfaces and any other sources of ignition. Do not smoke. Do not pierce or burn, even after use. Protect from sunlight. Do not expose to temperatures exceeding 50 ° C / 122 ° F

STERILIZATION METHODS

Device is supplied as non-sterile.

CLEANSING

This medical device is intended for repeated use, it is used until its contents are consumed.
After each use, wash the actuator in a warm soapy solution, then rinse in running water and wipe until completely dry.

OPERATIONAL, STORAGE AND TRANSPORTATION CONDITIONS

This medical device is transported by means of all types of transport in roofed vehicles in accordance with the requirements and rules for the transport of goods operating on vehicles of each type. Transportation conditions: At a temperature from 10 to 50 °C and a relative humidity of NMT 80%. During operation the device is resistant against influence of climatic aspects: temperatures from 10 to 50 °C and a relative humidity of NMT 80%. Store at the temperature from 10 to 50 °C.

SHELF-LIFE

3 years from the date of manufacture. Do not use the device after the expiration date indicated on the package.

ENVIRONMENTAL REQUIREMENTS

This medical device does not adversely affect humans or the environment when used, transported or stored.

LABELLING

Graphic symbols on the labelling of medical product

Graphic symbols	Description
	Manufacturer
	Read the package leaflet carefully before use!
	"Throw in the bin" - the device and/or its package are disposed with household waste at home.
	Complies with European Directive 93/42/EEC
	Batch number
	Expiry date
	Manufacturing date

PROCEDURE FOR DISPOSAL OR DESTRUCTION OF MEDICAL DEVICE

Special precautions for the destruction of an unused/used device are not required. The treatment and its

disposal must be carried out in accordance with the standards accepted in medical practice, as well as relevant local, state and federal laws and instructions. In medical institutions, device and/or its package are disposed of as Class A medical waste in accordance with SanPin 2.1.7.2790-10. At home, the device and/or its packaging are disposed of with household waste.

LIST OF INTERNATIONAL STANDARDS

Nasal spray MARIMER complies with the following standards:

Standard	Name
EN ISO 13485	Standard: Medical devices - QUALITY MANAGEMENT SYSTEMS - Requirements for regulatory purposes.
EN ISO 14971	Medical devices - Application of risk management to medical devices.
Council Directive 75/324/EEC of 20 May 1975	On the approximation of the laws of the Member States relating to aerosol dispensers.
EN ISO 10993-1	Biological evaluation of medical devices - Part 1 Evaluation and testing within a risk management process.
EN ISO 10993-5	Biological evaluation of medical devices - Part 5 Tests for in vitro cytotoxicity.
EN ISO 10993-10	Biological evaluation of medical devices - Part 10 Tests for irritation and skin sensitization.
EN 62366	Medical devices – Application of usability engineering to medical devices.
EN 1041	Information supplied by the manufacturer of medical device.
EN ISO 15223-1	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1 General requirements.

LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN

Does not contain materials of animal and (or) human origin.

INFORMATION ABOUT DRUG PRODUCTS, MEDICINAL PRODUCTS AND PHARMACEUTICAL SUBSTANCES

Does not contain drug products, medicinal products or pharmaceutical substances.

WARRANTY

This device is designed, manufactured, tested and packaged in compliance with all relevant requirements. The manufacturer guarantees the quality of the medical device until the expiration date, provided that its package is not damaged, and storage conditions are met.

RECLAMATION

Reclamations to be send to the address of the authorized representative of the manufacturer in the Russian Federation:

Limited liability company «Abbott Laboratories» («Abbott Laboratories» LLC)
125171, Moscow, Leningradskoe shosse, 16A, build. 1
Tel./fax: +7(495)-258-42-80/81

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