B|BRAUN

B. Braun Avitum AG
Division B. Braun Avitum

Schwarzenberger Weg 73-79
34212 Melsungen


To Federal Service on Surveillance in Healthcare
«Roszdravnadzor»

Place of issue: Melsungen


Date of issue: 2021-02-17

We, B. Braun Avitum AG, a company registered at: Schwarzenberger Weg 73 - 79, 34212 Melsungen, Germany, hereby certify authenticity of the attached document: «Instruction for use to the medical device «DiaStream® iQ» blood supply line kit for the Dialog iQ hemodialysis machine»

B. Braun Avitum AG

Dr. Holger Seeberg
Member of the Management Board B. Braun Avitum AG

Dr. Bettina Teschner
Senior Vice President Quality Management & RA

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
Carefully read these instructions for use before usage.

Refer also to the instructions for use of Dialog® iQ and of all medical devices/products that DiaStream® iQ bloodlines are used in combination with.

Product description

DiaStream® iQ comprise single-use extracorporeal A/V bloodlines. They represent one component of a hemodialysis treatment system.

The materials in direct or indirect contact with blood are: PP, PE, PVC, PI, ABS, TPE, EPTFE membrane.

All B. Braun A/V bloodlines do not contain natural rubber latex.

Color-coded components identify the arterial (red) and venous (blue) bloodlines.

Specific details on product configuration can be found on labelling.

Check symbols description in these instructions for use.

Refer to the product labeling for the volume of blood circuit.

Intended use

The DiaStream® iQ is intended to be used exclusively with the Dialog® iQ machine. DiaStream® iQ bloodlines in combination with Dialog® iQ machine are intended for hemodialysis treatments: Hemodialysis (HD), Hemodiafiltration (HDF), Hemofiltration (HF), sequential (SEQ)/ isolated ultrafiltration (ISO UF).

They are intended to provide the extracorporeal blood circuit from the patient's vascular access to the dialyzer and back in combination with the hemodialysis machines.

Indication for use

The DiaStream® iQ bloodlines are indicated for patients, who need hemodialysis treatments with chronic hemodialysis machines due to renal insufficiency in chronic kidney disease and/or acute kidney injury.

Contraindications

Therapy: Generally, contraindications for a hemodialysis treatments are known hypersensitivity to any material used and/or possibly the patient condition (clinical aspects, uncontrollable coagulation anomalies, etc.).

Product: Known hypersensitivity of the patient to any of the materials used in this product (see technical data).

Intended use environment

Hospitals, healthcare or limited-care facilities and intensive care units.

Emergency medication for the management of possible side effects of treatment must be readily available.

Patient population

Based on patient characteristics (cardiovascular status, hemodynamic stability, comorbidities, therapy tolerance, body size and weight, fluid and blood volume status, etc.) and clinical requirements the physician is responsible

- to assess carefully the risks caused by the amount of extracorporeal blood volume (especially low weight patients).
- to observe the weight limitations as indicated in the instructions for use of the hemodialysis machine.

Not for infants.

In pregnant women and nursing mothers, hemodialysis treatments should be used with particular caution. The responsible physician has to assess potential risk for mother and fetus in the prescription.

Intended users

Healthcare professionals trained and instructed for its appropriate use according to these instructions for use included in the training of Dialog® iQ machine.

In limited care settings patients/persons trained by healthcare professionals for operating the machine and all medical devices that the machine is used in combination with.

The responsible organization must ensure that the instructions for use are read and understood by all persons entrusted with any kind of work on or with the machine/product. The instructions for use must be permanently available to the user.

Side effects

Side effects related to the therapy – hemodialysis treatments

Refer to the instructions for use of Dialog® iQ machine.

Side effects related to the product – A/V bloodlines

In rare cases hypersensitivity reactions including anaphylactic (with possible shock and death) or anaphylactoid reactions to the A/V bloodlines, or other materials in the extracorporeal circuit, may occur during hemodialysis. Hypersensitivity reactions may cause mild to severe signs and symptoms including: malaise, nausea, headache, itching, flushing, hives, peripheral and facial swelling, erythema, ocular hyperemia, tingling sensation of the mouth and jaw, fever, leukopenia, hemolysis, anemia, hypotension, hypertension, tachycardia, arrhythmias, shortness of breath (dyspnea), wheezing, asthmatic reactions, bronchospasm, chest congestion, intradialytic pulmonary hypertension, reduced oxygen concentration and/or respiratory arrest, hemoconcentration, convulsion, unconsciousness, chronic low grade systemic inflammation, complement activation, immune dysregulation.

If severe hypersensitivity reactions occur, the dialysis must be discontinued and appropriate aggressive medical treatment for anaphylaxis must be initiated. The blood from the extracorporeal system must not be returned to the patient.

Patients with a history of hypersensitivity reactions or patients who have a history of being highly sensitive and allergic to a variety of substances should be carefully monitored during treatment.

Residual Risks

Residual risks in therapy – hemodialysis


Refer to the instructions for use of the Dialog® iQ machine.

Residual risks of the product – A/V bloodlines

- Blood loss to the environment may occur due to loose connection, disconnection or damage of Luer-Lock.
- Leakage of substitution fluid due to loose connection, disconnection or damage of Luer-Lock may result in incorrect patient weight control and hemoconcentration.
- The roller pump can provoke the release of micro-particles of PVC (Particle spallation) into the bloodstream, with possible long-term accumulation in lung, spleen and liver.

Warnings, cautions & precautions

- Handle the A/V bloodlines following aseptic technique.
- Remove protective connector caps only immediately before product connection to related devices used in the hemodialysis treatments
- Always connect B. Braun A/V bloodlines to vascular access devices with suitable Luer-Lock connectors. The manufacturer takes no responsibility for any use of products with incompatible features or layout. The use of other connectors differing from the above mentioned standard may cause problems such as stress cracks and leaks associated with the risk of air embolism and blood loss.
- Ensure the A/V line is correctly placed in the proper Safety Air Detector system of the respective dialysis machines, to prevent the risk of air infusion.
- Do not use needles larger than 21 gauge when accessing the injection ports.
- The injection port shall be disinfected prior to puncturing, using isopropyl alcohol 70% or chlorhexidine water based solution 2%. All the connectors shall be disinfected using only chlorhexidine water based solution 2%. The compatibility with any other disinfectant shall be determined prior to use.
- Administration of any drugs before or during treatment is the responsibility of the physician.
- Should the hydrophobic membrane of the pressure transducer protector (where present) come in contact with fluid, the A/V bloodline must be replaced.
- In case of occurrence of positive or negative pressure alarms during treatment, the POD membranes temporarily change position. In such event,
 1. do not manipulate the bloodline assembly on the machine
 2. do not disconnect the POD pressure monitoring lines for PA (negative pressure before blood pump) and PBE (positive pressure before dialyzer) from machine's PA or PBE pressure transducer
 3. check carefully for any cause of the alarm (e.g. arterial access occlusion, arterial line kinking, arterial clamp closed, ...)
 4. after having solved the alarm, the membranes will automatically return back in the correct position, but this is only possible if bloodline assembly on the machine has not been manipulated (see above)
 5. If the alarm was generated by a manipulation of the POD by the operator, the above-mentioned procedure will not be enough to restore the POD membrane position. In such case, proceed according to the following instructions:

Correct POD Membranes Positioning		
	Figure 1 PA / Arterial POD Correct membrane position	Figure 2 PBE POD Correct membrane position
Please ensure PODs are positioned vertically on the machine. Afterwards, please check the correct positioning of the membranes as displayed on the Figure 1 and Figure 2 on the right. The Figures 1 and 2, as well as the Figures 3 and 4, represents the POD profile sections. The thick black shape represents the membrane, the grey shape the blood.		
Incorrect POD Membranes Positioning		
Figure 3 PA / Arterial POD Incorrect membrane position	Figure 4 PBE POD Incorrect membrane position	
		
Scenarios of incorrect positioning of the PA / Arterial POD membrane:	Scenarios of incorrect positioning of the PBE POD membrane:	
1) PA / Arterial POD membrane is against the right wall with arterial pressure alarm on (refer to Figure 3); 2) No arterial pressure monitoring is shown on the dialysis machine.	1) PBE pressure POD membrane is against the left wall with PBE pressure alarm on (refer to figure 4); 2) No PBE pressure monitoring is shown on the dialysis machine.	
Instruction to rectify membrane positioning:	Instruction to rectify membrane positioning:	
1) Stop the blood pump; 2) Disconnect the monitoring line from arterial pressure sensor; 3) Keep the monitoring line disconnected until the membrane returns to the correct position (please refer to Figure 1); 4) Reconnect the monitoring line to the arterial pressure sensor; 5) Restart the blood pump. If PA alarm occurs, please check the patient arterial access. Should this not solve the alarm, please refer to the dialysis machine IFU.	1) Stop the blood pump; 2) Disconnect the monitoring line from the PBE pressure sensor; 3) Connect an empty sterile syringe with the plunger pulled back to the monitoring line. Push approximately 4 ml of air until the membrane returns to the correct position (please refer to Figure 2); 4) Clamp the monitoring line, disconnect the syringe and reconnect the monitoring line to the PBE pressure sensor. Thereafter, unclamp the monitoring line; 5) Restart the blood pump. If PBE alarm occurs, please refer to the dialysis machine IFU.	

Storage and handling conditions

Storage temperature between 0°C and 30°C.

Keep away from direct sunlight.

Do not use a sharp tool when opening the box or an individual package.

Keep A/V bloodlines stored in an undamaged packaging to ensure sterility. Do not use the device if packaging is damaged or protective caps are not in place. The expiry date only applies when the A/V bloodlines are correctly stored in an undamaged packaging. Do not use this product after the expiry date.

Sterility

Sterilized by irradiation. Sterile fluid path, non-pyrogenic fluid path.

Remove the A/V bloodlines from the sterile packaging only immediately before use.

Number of uses & duration of use

A/V bloodlines are intended for one dialysis treatment only. Re-use of single-use devices creates a potential risk to the patient or user. It may lead to contamination and/or infringed functionality, which may result in injury, illness or death of the user or the patient.

Do not use for more than 12 hours.

Clinical benefits:

Hemodialysis provides a life-saving treatment and includes the following clinical benefits:

- Removal of excess fluid
- Removal of uremic retention solutes and molecules
- Electrolyte removal and control
- Acid-Base adjustment

The DiaStream® iQ bloodlines contribute to all of the above being part of the extracorporeal blood circuit.

Performance characteristics

The A/V bloodlines are designed to withstand the maximum and minimum pressures and flow rates generated in use with the indicated dialysis machines.

NOTE: the actual blood flow rate might differ from the blood flow rate indicated by the machine. This difference might change with time and pressures.

Operating instructions

Set-up of the A/V bloodline systems

Please inspect the A/V bloodline set carefully for damage before usage.

a. Arterial set-up

- Inspect blood pump rotor before operation, looking for damages, deposits and dirt. To insert the blood pump segment, at first press upper part of multiconnector (green component) into blood pump housing, then lower part.
- Connect the arterial dialyzer connector to the dialyzer.

Pinching or any other damaging of the pump segment in the housing might cause pump segment rupture and consequent blood loss.

- Insert the arterial patient tubing into the HCT sensor and into the arterial clamp of the machine. If SNCO modality is performed, insert the arterial chamber into the respective holder of the machine.
- Attach the arterial patient tubing to the respective holder on the machine.
- Connect the monitoring line for the arterial pressure to the respective connector at the machine.
- Connect the monitoring line for the pre-filter pressure transducer (PBE) to the respective connector at the machine. If present, insert the PBE chamber into the respective holder.
- Connect the heparin line to the heparin syringe, if this anticoagulation method is performed. If heparin is not administered, clamp the heparin line to avoid air intake.

b. Venous set-up

- Insert the venous chamber into the respective holder.
- Ensure that the venous patient tube is properly inserted into the Safety Air Detector and into the venous clamp of the machine. Inaccurate insertion can lead to air embolism.

The Safety Air Detector will not detect air infused by a syringe through an access port downstream to the air detector.

- Attach the venous patient tubing to the respective holder on the machine.
- Connect the venous dialyzer connector to the dialyzer.
- Connect the monitoring line for venous pressure to the respective connector at the machine.

c. Additional instructions for HDF therapies

- Insert the pump segment of the substitution line into the respective pump rotor (refer to the handling instruction for blood pump rotor in Arterial set-up (a)).
- Connect lower tubing to the substitution port. Depending on post- or pre-dilution, connect the other tubing to the infusion port situated after the dialyzer (post-dilution) or to the infusion port situated after the blood pump and before the dialyzer (pre-dilution).

d. General set-up instructions

- To securely fix the transducer protectors (TP) of the pressure monitoring lines grab them by the wings and screw them firmly. Excessive tightening might cause breakage of the component. Do not grab them by the hydrophobic filter. Wet or leaking TP can lead to inaccurate pressure readings. In this case, TP must be replaced (if TP is removable), or bloodline (if TP is not removable).

Transducer protector (TP) / POD membrane ruptures might cause contamination of the machine with blood, leading to cross-contamination risk for patients. Refer to the instructions for use of the machine for necessary interventions by the technical service.

- Fix the tubing of the A/V bloodlines to the fixations of the machine to avoid their kinking and twisting during priming and treatment. Kinking, twisting and occlusion of the A/V bloodlines might result in hemolysis. Hemolysis may occur at elevated negative and positive pressures. Check the A/V bloodlines and patient's vascular access for potential kinking whenever an excessive pressure is observed.
 - Twisting and kinking of the pump segment/tubing set might cause improper occlusion, which may result in hemolysis and flow rate errors. Inspect the A/V bloodlines to detect visible defects and kinking.
 - Ensure that an appropriate dose of anti-coagulants are administered in order to prevent clotting. The indication, the selection and the dosage of anticoagulants is at the responsibility of the prescribing physician. Wrong concentration of anticoagulant may lead to injury (e.g. internal bleeding), illness or death of the patient, or may lead to thrombus formation in the capillaries of the dialyzer which may reduce its performance or result in clotting of the extracorporeal system.
 - To avoid air intake, keep the clamps of the infusion ports always closed during the treatment, except during infusion.
1. Priming of the blood circuit
- After setting-up make sure that all connections are tightly screwed.

Do not overtighten Luer-Lock connections, only tighten "finger-tight". Improper closure of Luer-Lock connectors, especially if too tight, can damage the connector leading to severe injury of the patient.

Too loose connection of Luer-Lock connectors can lead to patient harms, like air embolism or blood loss.

- Ensure that all clamps in the service lines are clamped, except the clamps on pressure monitoring lines and patient clamps.
 - Prime the extracorporeal system according to the instructions for use of the machine. The steps for loading and connection must be followed to ensure proper priming for efficient air and any remaining particle removal from the extracorporeal system.
 - Please comply with the dialyzer and dialysis machine manufacturer's instruction for use with regards to the required priming and rinsing volume. A minimum volume of 500ml is required to rinse the A/V Line in single path modality.
 - Check carefully for leaks during priming.
2. During treatment
- While connecting the bloodlines to the vascular access, make sure that the connection is tight.
 - Start and perform the treatment according to the instructions for use of the machine.
 - Retighten all connections at the beginning of the treatment, if necessary, ensuring to not overtighten the connections.
 - Check carefully and regularly for leaks during the treatment.
 - It is essential to regularly check the blood level in the bubble trap. The bubble trap must be filled to at least $\frac{3}{4}$.
 - Check regularly for clotting of the arterial pressure monitoring line and discontinue treatment if necessary.
 - During treatment make sure that the clamp of the heparin line, if present and used, is open.
 - Regularly check the heparin pump in order to verify the effective heparin infusion.
3. End of treatment
- Finish the treatment according to the instruction for use of the used machine.
 - Close all connections and dispose of the extracorporeal circuit as a closed system.

Disposal of the device

The used A/V bloodlines including all single use materials of the extracorporeal system has the potential to transmit infectious agents and shall be disposed appropriately immediately after use. The disposal must comply with the local laws and regulations.

Notice to the user

If any serious incident occurs in relation to the device this incident should be reported to B. Braun Avitum AG and to the responsible national competent authority.

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Обозначение детали в соответствии со схемой	Наименование детали
Артериальная кровопроводящая магистраль (Arterial Line) The diagram shows a complex assembly of tubes and connectors. It includes a main vertical tube on the left with various side branches and connectors. A coiled tube is shown in the upper right. Numerous numbered callouts (1-25) point to specific components like connectors, clamps, and tube sections. The diagram is a technical line drawing with some shading to indicate depth.	
1	Вставной закрытый колпачок (нейтральный) без ДЭГФ

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