
Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
www.medtronic.com
tel 001-763-514-4000
fax 001-763-514-4879


Declaration Letter

Medtronic Inc., a corporation organized under the laws of the United States, having its principal place of business at 710 Medtronic Parkway, Minneapolis, MN 55432, USA, hereby affirm that the attached documentation, IFU for MRI model 4196, is true and unaltered copies of the original document.

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The following list includes trademarks or registered trademarks of Medtronic in the United States and possibly in other countries. All other trademarks are the property of their respective owners.

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Attain, Attain Ability, Medtronic, SureScan

Contents

1	Description	3
2	Indications	4
3	Contraindications	4
4	Warnings and precautions	4
5	Potential adverse events	6
6	Implant procedure	7
7	Specifications	13
8	CE mark of conformity	15
9	Explanation of symbols on package labeling	15
10	Medtronic disclaimer of warranty	16
11	Service	16

1 Description

The Medtronic Attain Ability MRI SureScan 4196 steroid eluting, dual electrode, transvenous, over the wire, cardiac vein pacing lead is designed for pacing and sensing via a cardiac vein. The lead has been tested for use in the Magnetic Resonance Imaging (MRI) environment. All lead lengths for this lead model are MR Conditional. This lead contains 2 electrodes (designed to function as cathodes or negative terminals): the distal electrode, positioned at the distal tip, and the proximal electrode, positioned 21 mm from the tip. See Section 7.2, "Specifications drawing (nominal)", page 14 for a lead drawing.

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The Model 4196 lead can be positioned with the aid of a guide wire or with a stylet. If a stylet is used, use only the stylets packaged with the lead or in a stylet kit (downsized knob). Always use a stylet that is 3 cm shorter than the lead length listed on the IS-1 connector label.

Note: To implant the Model 4196 lead in a cardiac vein, a compatible delivery system is required. A compatible delivery system includes a guide catheter and either a hemostasis valve or an introducer valve that can be removed or that allows passage over the IS-1 connector. Contact a Medtronic representative for further information regarding compatible delivery systems.

1.1 Medtronic SureScan system

The Model 4196 lead is part of the Medtronic SureScan system. The SureScan system includes a Medtronic SureScan device connected to Medtronic SureScan leads. Labeling for SureScan system components displays the SureScan logo and the MR Conditional symbol.

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The MRI SureScan feature permits a mode of operation that allows a patient with a SureScan device to be safely scanned by an MRI machine while the device continues to provide appropriate pacing. When programmed to On, MRI SureScan operation disables arrhythmia detection, magnet mode, and all user-defined diagnostics. Before performing an MRI scan, refer to the SureScan MRI technical manual for important information about procedures and MRI-specific warnings and precautions.

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3 Contraindications

Before performing an MRI scan, refer to the SureScan MRI technical manual for MRI-specific contraindications. The Model 4196 lead is contraindicated as follows:

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Magnetic resonance imaging (MRI) – An MRI is a type of medical imaging that uses magnetic fields to create an internal view of the body. If certain criteria are met and the warnings and precautions provided by Medtronic are followed, patients with an MRI Conditional device and lead system are able to undergo an MRI scan; for details, refer to the SureScan MRI technical manual that Medtronic provides for an MRI Conditional device.

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- air embolism
- cardiac dissection
- cardiac perforation
- cardiac tamponade
- coronary sinus dissection
- death
- endocarditis
- erosion through the skin
- extracardiac muscle or nerve stimulation
- fibrillation or other arrhythmias
- heart block
- heart wall or vein wall rupture
- hematoma/seroma
- infection
- myocardial irritability

- myopotential sensing
- pericardial effusion
- pericardial rub
- pneumothorax
- body rejection phenomena including local tissue reaction
- threshold elevation
- thrombosis
- valve damage (particularly in fragile hearts)

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Implant techniques that may damage the lead include, but are not limited to, the following:

Implant techniques that may damage the lead	Possible effects on the lead	Corrective action to be considered
Forcing the lead through the introducer/delivery system	Electrode, conductor coil, and/or insulation damage	Replace the lead.
Use of too medial of an approach with venous introducer resulting in clavicle and first rib binding	Conductor coil fracture, insulation damage	Replace the lead.
Using too stiff a stylet	Conductor coil/insulation perforation	Replace the lead.
Puncturing the pericardium and/or tendon when using subclavian introducer approach	Conductor coil fracture, insulation damage	Replace the lead.
Advancing the lead through the veins without the stylet or guide wire fully inserted	Tip distortion and/or insulation perforation	Replace the lead.

Implant techniques that may damage the lead	Possible effects on the lead	Corrective action to be considered
Inserting the proximal end of the guide wire through the lead tip seal without using the guide wire insertion tool	Lead tip seal damage and/or conductor coil/insulation damage	Replace the lead.
Advancing a stylet tip beyond the distal end of the lead tip seal	Lead tip seal damage	Replace the lead.

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Note: To implant the Model 4196 lead in a cardiac vein, a compatible delivery system is required, such as a Medtronic delivery system. A compatible delivery system includes a guide catheter and a hemostasis valve or introducer valve which allows passage through or removal from an IS-1 connector. Contact a Medtronic representative for further information regarding compatible delivery systems.

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
Access the subclavian vein:

- Use a preferred method, based on professional experience, to access the subclavian vein.

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Figure 1.

1 Suggested entry site

- Introduce a J-shaped introducer guide wire and percutaneous introducer sheath.

6.4 Inserting the guide catheter assembly

Access the coronary sinus:

- Insert the guide catheter assembly.
- Advance the guide catheter to the right atrium.
- Achieve entry into the coronary sinus by rotating the guide catheter tip posteriorly and to the patient's left.

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Warning: Use only the stylets packaged with the lead or in a stylet kit (downsized knob). Always use a stylet that is 3 cm shorter than the lead length listed on the IS-1 connector label. Other stylets may extend beyond the lead tip causing lead tip seal damage or patient injury.

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
1. Insert a straight stylet into the lead to vary the shape of the distal end of the lead (Figure 2).

Note: When the stylet is fully inserted, the distal tip of the stylet does not reach the distal tip of the lead.

Figure 2.

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
The Model 4196 lead can be placed with the aid of a stylet or with a guide wire.

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

There are 2 techniques that may be used to advance the lead into a cardiac vein. The choice of technique is left to the discretion of the physician.

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Figure 6.


2. Use 1 or more of the following actions to direct the lead tip into the cardiac vein:

- Withdraw the stylet
- Advance the lead off of the stylet
- Remove the straight stylet. Manually shape the stylet into a curve. Insert the curved stylet into the lead. Rotate the lead body and stylet together.
- Rotate the curved tip of the lead into the vein.

3. Advance the lead into the cardiac vein (Figure 7).

Figure 7.

4. Advance the lead, with the stylet inserted, further into the cardiac vein. See Section 6.11, "Fixating the lead", page 11 for additional information.

Note: The proximal electrode should be placed in the cardiac vein, not in the coronary sinus.

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Figure 8.


2. Refer to 1 of the following procedures according to the type of stylet used.

- **Straight stylet:** If using a straight stylet, withdraw the stylet 1 to 2 cm and rotate the lead to align the curved lead tip with the cardiac vein ostium.
- **Curved stylet:** If using a curved stylet, withdraw the stylet 1 to 2 cm and rotate the stylet and lead body together to align the curved lead tip with the cardiac vein ostium.

3. Withdraw the lead and drag the curved lead tip over the ostium of the cardiac vein until the lead turns into the vein. Reinsert the stylet slightly to advance the lead into the cardiac vein (Figure 9).

Figure 9.

4. Advance the lead, with the stylet inserted, further into the cardiac vein. See Section 6.11, "Fixating the lead", page 11 for additional information.

Note: The proximal electrode should be placed in the selected vein, not the coronary sinus.

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Note: Medtronic recommends using guide wires 0.36 mm to 0.48 mm in diameter. Contact a Medtronic representative for further information about recommended guide wires.

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1. Select a guide wire. A more flexible guide wire is recommended if the patient has tortuous anatomy. If additional support is needed, use a firmer guide wire. **Note:** Guide wire firmness is determined by 2 factors: guide wire diameter and guide wire core design. A larger guide wire diameter may be firmer than a smaller guide wire diameter. However, guide wires with the same diameter may have different degrees of firmness. A "J" shaped guide wire or angled tip guide wire are recommended to aid in sub-selecting the cardiac vein and where more steerability is desired.

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

Figure 10.

- 1 Guide wire
- 2 Lead connector pin

3. Disengage the guide wire insertion tool from the lead connector pin.
4. Remove the guide wire insertion tool by sliding the tool off the end of the guide wire.
5. Position the guide wire steering handle (Figure 11) using the following steps :
 - a. Advance the guide wire steering handle over the proximal (rigid) end of the guide wire.
 - b. Tighten the guide wire steering handle onto the guide wire near the lead connector pin.

Figure 11.

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
Note: Grasp the lead and the guide wire insertion tool between your thumb and forefinger to prevent the lead from dislodging from the guide wire insertion tool while you are inserting the guide wire.

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There are 2 techniques that may be used to advance the lead into a cardiac vein using a guide wire. The choice of technique is left to the discretion of the physician.

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
1. Using venogram images as reference, advance the lead into the cardiac vein (Figure 15). The curved distal tip of the lead


may be used to aid in subselecting cardiac veins within the coronary sinus by slightly retracting the guide wire inside the lead lumen.

Figure 15.

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
Figure 17.


Note: If the cardiac vein is large, it may be necessary to position the lead in a smaller cardiac vein in order to achieve lead tip fixation. Both electrodes must touch the tissue.

6.12 Taking electrical measurements

The Model 4196 lead was designed to provide unipolar pacing via selectable cathodes. Recommended pacing vectors are shown in Figure 18.

Figure 18. Recommended Model 4196 lead pacing vectors

Warning: This lead has not been designed for pacing in the bipolar configuration. It utilizes 2 equal-sized electrodes. Use of this lead in a bipolar configuration may result in higher pacing thresholds or anodal stimulation.

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1. Attach a surgical cable to the lead connector pin or ring depending on which electrode you want to use for measurements.
Note: Testing the lead in a unipolar configuration requires the use of an indifferent electrode.
2. Use an implant support instrument, such as a pacing system analyzer, for obtaining electrical measurements. For information on the use of the implant support instrument, consult the product documentation for that device.
Satisfactory lead placement is indicated by low stimulation thresholds and adequate sensing of intracardiac signal amplitudes.
 - A low stimulation threshold provides a desirable safety margin, allowing for a possible rise in thresholds that may occur within 2 months following implant.

- Adequate sensing amplitudes ensure that the lead is properly sensing intrinsic cardiac signals. Minimum signal requirements depend on the sensitivity capabilities of the device. Acceptable acute signal amplitudes for the lead must be greater than the minimum device sensing capabilities. Be sure to include an adequate safety margin to account for lead maturity.

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1. If used, remove the guide wire and guide wire insertion tool. Replace the guide wire with a straight stylet (downsized knob). Insert the straight stylet into the lead to the mid-coronary sinus.
2. Remove the guide catheter from the lead. See the delivery system product documentation for details.
Note: For Medtronic slidable delivery systems, use a slitter compatible with a 1.3 mm (4 French) lead body.
3. Carefully and completely remove the stylet. When removing the stylet, grip the lead firmly just distal of the connector pin; this will help prevent possible lead tip dislodgement.
4. Repeat the electrical measurements. See Section 6.12, "Taking electrical measurements", page 12.

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- Use only nonabsorbable sutures to anchor the lead.
- Do not attempt to remove or cut the anchoring sleeve.
- Do not use the anchoring sleeve tabs for suturing.
- During lead anchoring, take care to avoid dislodging the lead tip.
- Do not secure sutures so tightly that they damage the vein, lead, or anchoring sleeve (Figure 19).
- Do not tie a suture directly to the lead body (Figure 19).

Figure 19.

Anchor the lead using all 3 grooves:

1. Position the anchoring sleeve against or near the vein.
2. Secure the anchoring sleeve to the lead body by tying a suture firmly in each of the 3 grooves (Figure 20).

Figure 20.

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
- Ensure that the leads do not leave the device at an acute angle.
- Do not grip the lead or device with surgical instruments.
- Do not coil the lead. Coiling the lead can twist the lead body and may result in lead dislodgement (Figure 21).

Figure 21.

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1. To prevent undesirable twisting of the lead body, rotate the device to loosely wrap the excess lead length (Figure 22).

Figure 22.

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7 Specifications

7.1 Specifications (nominal)

Parameter	Model 4198
Type	Dual electrode
Chamber paced	Left ventricle
Length	78 and 88 cm
Connector	IS-1 BI
Material	Conductor: Si-Pt coated 25% Ag-core-MP35N
	Insulators: Polyurethane (outer) Si-polyimide (Si-Pi) (inner) ^a
	Ring electrode: Platinum iridium alloy with titanium nitride coating
	Tip electrode: Platinum iridium alloy with titanium nitride coating
	Connector pin: Stainless steel
	Connector ring: Stainless steel
	Molded tip seal: Silicone rubber (with steroid)


English 13

Parameter	Model 4196
Tip electrode configuration	Tapered, annular, titanium-nitride-coated, steroid-eluting
Ring electrode configuration	Radiused, titanium-nitride-coated, steroid-eluting
Diameter	Lead body: 1.3 mm (4 French)
	Tip electrode: 1.53 mm (4.6 French)
	Ring electrode: 1.70 mm (5.1 French)
Medtronic delivery system (recommended inner diameter)	1.90 mm (5.7 French)
Diagnostic guide wire (recommended diameter)	0.36 mm to 0.46 mm
Electrode surface area	Tip: 5.8 mm ²
	Ring: 5.8 mm ²
Conductor resistance	Tip: 38 ± 10 Ω (78 cm)
	43 ± 10 Ω (88 cm)
	Ring: 35 ± 10 Ω (78 cm)
	40 ± 10 Ω (88 cm)
Distance between electrodes	21 mm
Steroid	Dexamethasone acetate
Combined total of steroid	<1.0 mg

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8 CE mark of conformity

The following CE mark of conformity applies to the Model 4196 lead.

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Symbol	Explanation
	Conformité Européenne (European Conformity). This symbol means that the device fully complies with applicable European Union acts.
	Do not use if package is damaged
	Do not reuse
	Upper limit of temperature
	Sterilized using ethylene oxide
	Consult instructions for use
	Consult instructions for use at this website
	Date of manufacture
	Manufacturer
	Use by
	Reorder number
	Authorized representative in the European Community
	Serial number

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
53010

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Declaration Letter


Medtronic Inc., a corporation organized under the laws of the United States, having its principal place of business at 710 Medtronic Parkway, Minneapolis, MN 55432, USA, hereby or affirm that the attached documentation, IFU for MRI model 4296, is true and unaltered copies of the original document.

Name: Shellee Lamb
Title: Sr. Regulatory Affairs Specialist

Oct 4, 2019

Date

Chidi A. Bigelbach
10/4/19

Medtronic

ATTAIN ABILITY™ PLUS MRI SURESCAN™ 4296

Steroid eluting, dual electrode, transvenous, over the wire, cardiac vein pacing lead

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Attain, Attain Ability, Attain Hybrid, Medtronic, SureScan

Contents

1	Description	3
2	Indications	4
3	Contraindications	4
4	Warnings and precautions	4
5	Potential adverse events	6
6	Implant procedure	7
7	Specifications	14
8	CE mark of conformity	15
9	Explanation of symbols on package labeling	15
10	Medtronic disclaimer of warranty	16
11	Service	16

1 Description

The Medtronic Attain Ability Plus MRI SureScan 4296 steroid eluting, dual electrode, transvenous, over the wire, cardiac vein pacing lead is designed for pacing and sensing via a cardiac vein. The lead has been tested for use in the Magnetic Resonance Imaging (MRI) environment. All lead lengths for this lead model are MR Conditional. This lead contains 2 electrodes (designed to function as cathodes or negative terminals): the distal electrode, positioned at the distal tip, and the proximal electrode, positioned 21 mm from the tip. See Section 7.2, "Specifications drawing (nominal)", page 15 for a lead drawing.

The distal tip electrode allows a guide wire to pass through to aid in cardiac vein selection. The distal tip electrode also contains a dual-function steroid-loaded silicone rubber membrane MCRD¹. This membrane seals the lead inner lumen to reduce blood ingress, and the MCRD elutes steroid to reduce inflammatory response within the cardiac vein. The proximal electrode is positioned between the 2 distal curves of the lead to facilitate contact with the cardiac veins. It also contains an MCRD for elution of steroid. The outer insulation of the lead is polyurethane and the inner insulation is SI-polyimide (SI-PI)². The SI-PI is applied as a coating to the conductor wire prior to coiling. The IS-1³ connector pin has a drilled inner lumen to improve guide wire passage.

The MCRDs contain a combined-total dosage of <1.0 mg of dexamethasone acetate steroid. Upon exposure to body fluids, the steroid elutes from the MCRDs. The steroid suppresses the inflammatory response that is believed to cause threshold rise typically associated with implanted pacing electrodes.

The Model 4296 lead can be positioned with the aid of a guide wire or with a stylet. If a stylet is used, use only the stylets packaged with the lead or in a stylet kit (downsized knob). Always use a stylet that is 3 cm shorter than the lead length listed on the IS-1 connector label.

Note: To implant the Model 4296 lead in a cardiac vein, a compatible delivery system is required. A compatible delivery system includes a guide catheter and either a hemostasis valve or an introducer valve that can be removed or that allows passage over the IS-1 connector. Contact a Medtronic representative for further information regarding compatible delivery systems.

1.1 Medtronic SureScan system

The Model 4296 lead is part of the Medtronic SureScan system. The SureScan system includes a Medtronic SureScan device connected to Medtronic SureScan leads. Labeling for SureScan system components displays the SureScan logo and the MR Conditional symbol.

SureScan logo

MR Conditional symbol. The Medtronic SureScan system is MR Conditional and is designed to allow implanted patients to undergo an MRI scan under the specified MRI conditions for use.

The MRI SureScan feature permits a mode of operation that allows a patient with a SureScan device to be safely scanned by an MRI machine while the device continues to provide appropriate pacing. When programmed to On, MRI SureScan operation disables arrhythmia detection, magnet mode, and all user-defined diagnostics. **Before performing an MRI scan, refer to the SureScan MRI technical manual for important information about procedures and MRI-specific warnings and precautions.**

1.2 Contents of package

Leads and accessories are supplied sterile. Each package contains the following items:

- 1 lead with anchoring sleeve
- 1 guide wire insertion tool
- 1 guide wire clip
- 1 guide wire steering handle
- extra stylets
- product documentation

¹ Monolithic controlled release device.

² Technology developed by NASA.

³ IS-1 refers to an International Connector Standard (ISO 5841-3) whereby pulse generators and leads so designated are assured of a basic mechanical fit.

English 3

1.3 Accessory descriptions

Dispose of all single-use accessories according to local environmental requirements.

Anchoring sleeve – An anchoring sleeve secures the lead to prevent it from moving and protects the lead insulation and conductors from damage caused by tight sutures.

Guide wire clip – A guide wire clip secures the excess guide wire and helps to protect and maintain the sterility of the guide wire.

Guide wire insertion tool – A guide wire insertion tool provides additional control when inserting a guide wire into the lead connector pin or the lead tip.

Guide wire steering handle – A guide wire steering handle is used only with guide wires 0.46 mm or less in diameter. The steering handle provides additional control and steerability of the guide wire.

Stylet – A stylet provides additional stiffness and controlled flexibility for maneuvering the lead into position. Each stylet knob is labeled with the stylet diameter and corresponding lead length.

2 Indications

The lead has application as part of a Medtronic biventricular pacing system.

3 Contraindications

Before performing an MRI scan, refer to the SureScan MRI technical manual for MRI-specific contraindications.

The Model 4296 lead is contraindicated as follows:

Coronary vasculature – This lead is contraindicated for patients with coronary venous vasculature that is inadequate for lead placement, as indicated by venogram.

Steroid use – The lead is contraindicated in patients for whom a single dose of 1.0 mg of dexamethasone acetate may be contraindicated.

4 Warnings and precautions

A complete SureScan system is required for use in the MRI environment. Before performing an MRI scan, refer to the SureScan MRI technical manual for MRI-specific warnings and precautions.

Note: Medical procedure warnings and precautions that pertain to the Medtronic implanted system are provided in the manual that is packaged with the device or on the Medtronic Manual Library website (www.medtronic.com/manuals).

Inspecting the sterile package – Inspect the sterile package with care before opening it.

- Contact a Medtronic representative if the seal or package is damaged.
- Do not store this product above 40 °C.
- Do not use the product after its expiration date.

Sterilization – Medtronic has sterilized the package contents with ethylene oxide before shipment. This lead is for single use only and is not intended to be resterilized.

Single use – The lead and accessories are for single use only.

Necessary hospital equipment – Keep external defibrillation equipment nearby for immediate use during acute lead system testing, the implant procedure, or whenever arrhythmias are possible or intentionally induced during post-implant testing.

Line-powered and battery-powered equipment – An implanted lead forms a direct current path to the myocardium. During lead implant and testing, use only battery-powered equipment or line-powered equipment specifically designed for this purpose to protect against fibrillation that may be caused by alternating currents. Line-powered equipment used in the vicinity of the patient must be properly grounded. Lead connector pins must be insulated from any leakage currents that may arise from line-powered equipment.

Concurrent devices – Output pulses, especially from unipolar devices, may adversely affect device sensing capabilities. If a patient requires a separate stimulation device, either permanent or temporary, allow enough space between the leads of the separate systems to avoid interference in the sensing capabilities of the devices. Previously implanted pulse generators and implantable cardioverter defibrillators should generally be explanted.

Steroid use – It has not been determined whether the warnings, precautions, or complications usually associated with injectable dexamethasone acetate apply to the use of this highly localized, controlled-release lead.

Pregnancy – Dexamethasone acetate has been shown to be teratogenic in many species when given in doses equivalent to the human dose. There are no adequate and well-controlled studies in pregnant women. Dexamethasone acetate should be used during pregnancy only if the potential benefit justifies the potential risk to the fetus. Studies in mice, rats, and rabbits have shown that adrenocorticoids increase the incidence of cleft palate, placental insufficiency, and spontaneous abortions, and can decrease the intrauterine growth rate.

Nursing mothers – Systemically administered corticosteroids appear in human milk and could suppress growth, interfere with endogenous corticosteroid production, or cause other untoward effects in nursing infants. Because of the potential for serious adverse reactions in nursing infants from corticosteroids, a decision should be made whether to discontinue nursing or to use a non-steroidal lead, taking into account the importance of the lead and the drug to the mother.

Handling steroid MCRDs⁴ – Avoid reducing the amount of steroid available before implanting the lead. Reducing the

available amount of steroid may adversely affect low-threshold performance.

- Do not allow the electrode surfaces to come in contact with surface contaminants.
- Do not wipe or immerse the electrodes in fluid, except blood, at the time of implant.

Bipolar pacing – This lead has not been designed for pacing in the bipolar configuration. It utilizes 2 equal-sized electrodes. Use of this lead in a bipolar configuration may result in higher pacing thresholds or anodal stimulation.

Handling the stylet – Handle the stylet with care at all times.

- To minimize the likelihood of trauma to the vein and to maintain lead flexibility while advancing the lead through the vein, keep the stylet withdrawn 1 to 2 cm or select a more flexible stylet.
- Do not use excessive force or surgical instruments when inserting a stylet.
- Avoid overbending, kinking, or blood contact on stylets.
- Use a new stylet when blood or other fluids accumulate on the stylet. Accumulated fluids may cause damage to the lead or difficulty in passing the stylet through the lead.
- Curving the distal end of the stylet prior to insertion into the lead will achieve a curvature at the distal end of the lead. Do not use a sharp object to impart a curve to the distal end of the stylet.

Handling the guide wire – Handle the guide wire with care at all times.

- Do not insert the proximal end of the guide wire through the lead tip seal without using the guide wire insertion tool. Inserting the guide wire without the guide wire insertion tool could cause damage to the lead tip seal or to the conductor core or insulation.
- Damage to a guide wire may prevent the guide wire from performing with accurate torque response and control and may cause vessel damage. For additional information about vessel damage and other potential adverse events, refer to the technical manual packaged with the appropriate guide wire.
- If the distal end of the guide wire becomes severely kinked or twisted, it may be difficult to withdraw it back through the lead. Therefore, if there is an indication that the distal end of the guide wire has become damaged, or if there is significant resistance in guide wire passage, remove the lead and guide wire together as a unit. Remove the guide wire from the lead and insert a new guide wire into the lead. Do not use excessive force to retract the guide wire from the lead. Refer to the product documentation packaged with the guide wire for additional information.

Handling the lead – Handle the lead with care at all times:

- Rust stylets are not recommended with this lead due to the risk of conductor coil or insulation perforation.
- If a stylet is used for lead positioning, use only the stylets packaged with the lead or in a stylet kit (downsized knob).

Other stylets may extend beyond the lead tip causing lead tip seal damage or patient injury.

- If the lead is damaged, do not implant it. Return the lead to a Medtronic representative.
- Protect the lead from materials that shed particles such as lint and dust. Lead insulators attract these particles.
- Handle the lead with sterile surgical gloves that have been rinsed in sterile water or a comparable substance.
- Do not severely bend, kink, or stretch the lead.
- Do not use surgical instruments to grasp the lead or connector pin.
- Do not immerse leads in mineral oil, silicone oil, or any other liquid, except blood, at the time of implant.
- Use an anchoring sleeve with all leads. Ensure that the anchoring sleeve is positioned close to the lead connector pin, to prevent inadvertent passage of the sleeve into the vein. If wiping the lead is necessary prior to insertion, ensure that the anchoring sleeve remains in position.
- Do not force the guide catheter or leads if significant resistance is encountered. Use of guide catheters or leads may cause trauma to the heart.

Chronic repositioning or removal – Chronic repositioning or removal of leads may be difficult because of fibrotic tissue development. The clinical study was not designed to evaluate the removal of left ventricular leads from the coronary venous vasculature. If a lead must be removed or repositioned, proceed with extreme caution. Return all removed leads to Medtronic.

- Rust stylets are not recommended with this lead due to the risk of conductor core/insulation perforation.
- Verify lead length on the IS-1 label on the connector to choose an appropriate stylet kit (downsized knob) length when repositioning. Always choose a stylet kit (downsized knob) 3 cm shorter than the lead length. For example, choose a stylet kit (downsized knob) with stylets 75 cm long for a lead 78 cm long.
- Lead removal may result in avulsion of the endocardium, valve, or vein.
- Lead junctions may separate, leaving the lead tip and bare wire in the heart or vein.
- Chronic repositioning may adversely affect the low-threshold performance of a steroid eluting lead.
- Cap abandoned leads to avoid transmitting electrical signals.
- For leads that have been severed, seal the remaining lead end and suture the lead to adjacent tissue.
- If a lead is removed and repositioned, inspect it carefully for insulator or conductor coil damage prior to repositioning.

Magnetic resonance imaging (MRI) – An MRI is a type of medical imaging that uses magnetic fields to create an internal view of the body. If certain criteria are met and the warnings and precautions provided by Medtronic are followed, patients with an MR Conditional device and lead system are able to undergo an MRI scan; for details, refer to the SureScan MRI technical manual that Medtronic provides for an MR Conditional device.

⁴ Monolithic controlled release device

Diathermy treatment (including therapeutic ultrasound) –

Diathermy is a treatment that involves the therapeutic heating of body tissues. Diathermy treatments include high frequency, short wave, microwave, and therapeutic ultrasound. Except for therapeutic ultrasound, do not use diathermy treatments on cardiac device patients. Diathermy treatments may result in serious injury or damage to an implanted device and lead system. Therapeutic ultrasound (including physiotherapy, high intensity therapeutic ultrasound, and high intensity focused ultrasound), is the use of ultrasound at higher energies than diagnostic ultrasound to bring heat or agitation into the body. Therapeutic ultrasound is acceptable if treatment is performed with a minimum separation distance of 15 cm between the applicator and the implanted device and lead system, as long as the ultrasonic beam is pointing away from the device and lead system.

External defibrillation and cardioversion – External defibrillation and cardioversion are therapies that deliver an electrical shock to the heart to convert an abnormal heart rhythm to a normal rhythm.

Medtronic cardiac devices are designed to withstand exposure to external defibrillation and cardioversion. While damage to an implanted system from an external shock is rare, the probability increases with increased energy levels. These procedures may also temporarily or permanently elevate pacing thresholds or temporarily or permanently damage the myocardium. If external defibrillation or cardioversion are required, consider the following precautions:

- Use the lowest clinically appropriate energy.
- Position the patches or paddles a minimum of 15 cm away from the device.
- Position the patches or paddles perpendicular to the device and lead system.
- If an external defibrillation or cardioversion is delivered within 15 cm of the device, use a Medtronic programmer to evaluate the device and lead system.

5 Potential adverse events

The potential adverse events related to the use of transvenous leads include, but are not limited to, the following patient-related conditions:

- air embolism
- cardiac dissection
- cardiac perforation
- cardiac tamponade
- coronary sinus dissection
- death
- endocarditis and pericarditis
- erosion through the skin
- extracardiac muscle or nerve stimulation
- fibrillation or other arrhythmias
- heart block
- heart wall or vein wall rupture
- hematoma/seroma
- infection
- myocardial irritability

- myopotential sensing
- pericardial effusion
- pericardial rub
- pneumothorax
- body rejection phenomena including local tissue reaction
- threshold elevation
- thrombosis
- thrombotic embolism
- valve damage (particularly in fragile hearts)

Additional potential adverse events related to the lead and the programmed parameters include, but are not limited to, the following:

Potential adverse event	Indicator of potential adverse event	Corrective action to be considered
Lead dislodgement ^a	Intermittent or continuous loss of capture or sensing ^a	Reposition the lead.
Lead dislodgement ^a	Intermittent or continuous oversensing	Reposition the lead.
Lead conductor fracture or insulation failure	Intermittent or continuous loss of capture or sensing ^a	Replace the lead.
Threshold elevation or exit block	Loss of capture ^a	Adjust the implantable device output. Replace or reposition the lead.

^a Transient loss of capture or sensing may occur following surgery until lead stabilization takes place. If stabilization does not occur, lead dislodgement may be suspected.

Implant techniques that may damage the lead include, but are not limited to, the following techniques:

Implant techniques that may damage the lead	Possible effects on the lead	Corrective action to be considered
Forcing the lead through the introducer/delivery system	Electrode, conductor coil, and/or insulation damage	Replace the lead.
Use of too medial of an approach with venous introducer resulting in clavicle and first rib binding	Conductor coil fracture, insulation damage	Replace the lead.
Using too stiff a stylet	Conductor coil/insulation perforation	Replace the lead.
Puncturing the pericardium and/or tendon when using subclavian introducer approach	Conductor coil fracture, insulation damage	Replace the lead.
Advancing the lead through the non-coronary central access veins without the stylet or guide wire fully inserted	Tip distortion and/or insulation perforation	Replace the lead.

6 English

Implant techniques that may damage the lead

Inserting the proximal end of the guide wire through the lead tip seal without using the guide wire insertion tool

Advancing a stylet tip beyond the distal end of the lead tip seal

Possible effects on the lead

Lead tip seal damage and/or conductor coil/insulation damage

Lead tip seal damage

Corrective action to be considered

Replace the lead.

Replace the lead.

6 Implant procedure

Warning: Before implanting a SureScan system, consider the risks associated with removing previously implanted leads. Abandoned leads or previously implanted leads not tested for MRI compatibility compromise the ability to safely scan the SureScan system during MRI scans.

Warning: Do not force the guide catheter or lead if significant resistance is encountered. Use of guide catheters or leads may cause trauma to the heart.

Note: To implant the Model 4296 lead in a cardiac vein, a compatible delivery system is required, such as a Medtronic delivery system. A compatible delivery system includes a guide catheter and a hemostasis valve or introducer valve which allows passage through or removal from an IS-1 connector. Contact a Medtronic representative for further information regarding compatible delivery systems.

Proper surgical procedures and sterile techniques are the responsibility of the medical professional. The implant procedures described in this manual are provided for information only. Each physician must apply the information in these instructions according to professional medical training and experience.

6.1 Placing the right ventricular lead

Note: When deciding which ventricular lead to place first, consider the ease of coronary sinus cannulation and the need for backup pacing.

- When a left ventricular lead is implanted first:
 - It may be easier to cannulate the coronary sinus without other leads already implanted.
 - Additional hospital equipment may be necessary to provide backup pacing.
- When a right ventricular lead is implanted first:
 - A right ventricular lead may be used to provide backup pacing.
 - It may be more difficult to cannulate the coronary sinus with a right ventricular lead already implanted.

6.2 Preparing the delivery system

Prepare the delivery system for lead implantation according to the instructions in the product documentation packaged with the delivery system.

6.3 Accessing the subclavian vein

Warning: Backup pacing should be readily available during implant. Use of the delivery system or leads may cause heart block.


Access the subclavian vein:

1. Use a preferred method, based on professional experience, to access the subclavian vein.

Caution: Certain anatomical abnormalities, such as thoracic outlet syndrome, may precipitate pinching and subsequent fracture of the lead.

Caution: Insertion should be done as far lateral as possible to avoid clamping the lead body between the clavicle and the first rib (Figure 1).

Figure 1.

1 Suggested entry site

2. Introduce a J-shaped introducer guide wire and percutaneous introducer sheath.

6.4 Inserting the guide catheter assembly

Access the coronary sinus:

1. Insert the guide catheter assembly.
2. Advance the guide catheter to the right atrium.
3. Achieve entry into the coronary sinus by rotating the guide catheter tip posteriorly and to the patient's left.

See the delivery system product documentation for additional information.

6.5 Obtaining venograms

Before placing a lead in the coronary sinus, obtain venograms. Venograms are recommended to assess a route for passage and a site for final placement based on the size, shape, location, and tortuosity of the veins. Also, venograms may be useful in identifying suspected coronary sinus trauma. For information on obtaining a venogram by using a venogram balloon catheter, see the product documentation packaged with an appropriate venogram balloon catheter.

English 7

6.6 Inserting the lead into the delivery system

Warning: If a stylet is used for lead positioning, use only the stylets packaged with the lead or in a stylet kit (downsized knob). Always use a stylet that is 3 cm shorter than the lead length listed on the IS-1 connector label. Other stylets may extend beyond the lead tip causing lead tip seal damage or injury or perforation of the cardiac vein or heart.

Warning: Do not force the lead if significant resistance is encountered during lead passage. The use of guide catheters or leads may cause trauma to the heart.


Caution: Use care when handling the lead during insertion.

- Do not severely bend, kink, or stretch the lead.
- Do not use surgical instruments to grasp the lead or connector pin.

Use the following steps to insert the lead:

1. Insert a straight stylet into the lead to vary the shape of the distal end of the lead (Figure 2).
Note: When the stylet is fully inserted, the distal tip of the stylet does not reach the distal tip of the lead.

Figure 2.

- 1 Stylet fully withdrawn
- 2 Stylet fully inserted

2. Insert the lead into the delivery system. See the delivery system product documentation for additional information.


6.7 Selecting the lead placement method

The Model 4296 lead can be placed with the aid of a stylet or with a guide wire.

To determine the placement method, review the venogram to determine if stylet or guide wire delivery should be used.

Stylet delivery – If the patient's anatomy features gentle vein angulation off of the coronary sinus and the cardiac vein branch is not tortuous (Figure 3), use a stylet for lead delivery.


Figure 3.

- 1 Coronary sinus
- 2 Cardiac vein

Guide wire delivery – If the patient's anatomy features acute vein angulation off of the coronary sinus and the cardiac vein branch is tortuous (Figure 4), use a guide wire for lead delivery.

Figure 4.

- 1 Tortuous cardiac vein branch with gentle angulation from the coronary sinus
- 2 Tortuous cardiac vein branch with acute angulation from the coronary sinus

6.8 Placing the lead using a stylet

Warning: Do not force the lead if significant resistance is encountered during lead passage.

Warning: To minimize the likelihood of trauma to the vein and to maintain lead flexibility while advancing the lead through the vein, keep the stylet withdrawn 1 to 2 cm or select a more flexible stylet.

Caution: To avoid damage to the stylet, do not use a sharp object to impart a curve to the distal end of the stylet (Figure 5).

Figure 5.

Note: If it is difficult to advance the stylet around a bend, consider changing the stylet. More flexible stylets are recommended for tortuous anatomies. Firmer stylets are recommended where additional support is needed.

There are many techniques that may be used to advance the lead into a cardiac vein using a stylet. The choice of technique is left to the discretion of the physician. Two possible techniques are described below.

6.8.1 First technique to advance the lead into a cardiac vein using a stylet

Use the following steps to advance the lead into a cardiac vein:


1. Advance the lead up to, but not past, the cardiac vein ostium (Figure 6).

Figure 6.

2. Use 1 or more of the following actions to direct the lead tip into the cardiac vein:
 - Withdraw the stylet
 - Advance the lead off of the stylet
 - Remove the straight stylet. Manually shape the stylet into a curve. Insert the curved stylet into the lead. Rotate the lead body and stylet together.
 - Rotate the curved tip of the lead into the vein.
3. Advance the lead into the cardiac vein (Figure 7).

Figure 7.

4. Advance the lead, with the stylet inserted, further into the cardiac vein. See Section 6.11, "Fixating the lead", page 11 for additional information.
- Note:** The proximal electrode should be placed in the cardiac vein, not in the coronary sinus.

6.8.2 Second technique to advance the lead into a cardiac vein using a stylet

Use the following steps to advance the lead into a cardiac vein:


1. Advance the lead past the ostium of the cardiac vein (Figure 8).

Figure 8.

2. Refer to 1 of the following procedures according to the type of stylet used.
 - **Straight stylet:** If using a straight stylet, withdraw the stylet 1 to 2 cm and rotate the lead to align the curved lead tip with the cardiac vein ostium.
 - **Curved stylet:** If using a curved stylet, withdraw the stylet 1 to 2 cm and rotate the stylet and lead body together to align the curved lead tip with the cardiac vein ostium.
3. Withdraw the lead and drag the curved lead tip over the ostium of the cardiac vein until the lead turns into the vein. Reinsert the stylet slightly to advance the lead into the selected cardiac vein (Figure 9).

Figure 9.

4. Advance the lead, with the stylet inserted, further into the cardiac vein. See Section 6.11, "Fixating the lead", page 11 for additional information.
- Note:** The proximal electrode should be placed in the selected vein, not the coronary sinus.

6.9 Preparing the guide wire

Warning: Do not insert the proximal end of the guide wire through the lead tip seal without using the guide wire insertion tool. Inserting the guide wire without the guide wire insertion tool could cause damage to the lead tip seal or to the conductor core or insulation.

Warning: Damage to a guide wire may prevent the guide wire from performing with accurate torque response and control and may cause vessel damage. For additional information about vessel damage and other potential adverse events, refer to the technical manual packaged with the appropriate guide wire.

Caution: Use care when positioning the guide wire. Refer to the product documentation packaged with the guide wire for additional information.

Note: Medtronic recommends using guide wires 0.36 mm to 0.46 mm in diameter. Contact a Medtronic representative for further information about recommended guide wires.

Use the following steps to prepare the guide wire for use:

1. Select a guide wire. A more flexible guide wire is recommended if the patient has tortuous anatomy. If additional support is needed, use a firmer guide wire.
Note: Guide wire firmness is determined by 2 factors: guide wire diameter and guide wire core design. A larger guide wire diameter may be firmer than a smaller guide wire diameter. However, guide wires with the same diameter may have different degrees of firmness. A "J" shaped guide wire or angled tip guide wire are recommended to aid in sub-selecting the cardiac vein and where more steerability is desired.

Note: Consider soaking the guide wire in a heparin solution prior to insertion to minimize the risk of thrombus formation during use.

Note: The insertion tool steps below, Step 2 through Step 4, do not apply when using the Attain Hybrid guide wire. The tool cannot be removed over the wire, due to the attached proximal knob. When following Step 5 and Step 6, the torque tool must be preloaded onto the Attain Hybrid guide wire, as the proximal attached knob prohibits backloading of the wire into the steering handle. The preloading step with the torque tool must be completed before loading the guide wire into the connector pin on the proximal end of the lead.


2. Insert the guide wire into the lead by placing the distal (flexible) end of the guide wire into the lead connector pin using the guide wire insertion tool included in the package (Figure 10).

Note: To prevent the lead from dislodging from the guide wire insertion tool, grasp the lead and guide wire insertion tool firmly between thumb and forefinger.

Caution: To minimize the risk of damaging the guide wire, be sure the flexible section of the guide wire is fully inserted into the lead before removing the guide wire insertion tool from the lead.

Note: Be sure to remove the guide wire insertion tool prior to lead implant.

Figure 10.


- 1 Guide wire
- 2 Lead connector pin

3. Disengage the guide wire insertion tool from the lead connector pin.
4. Remove the guide wire insertion tool by sliding the tool off the end of the guide wire.

5. Position the guide wire steering handle (Figure 11) using the following steps:

- a. Advance the guide wire steering handle over the proximal (rigid) end of the guide wire.
- b. Tighten the guide wire steering handle onto the guide wire near the lead connector pin.

Figure 11.

- 1 Guide wire steering handle
- 2 Lead connector pin

6. Attach the guide wire clip to the guide wire and secure it within the sterile field. Medtronic recommends securing the guide wire clip to the patient's sterile surgical drape.

As an alternate approach in situations where the guide wire is already in place, the lead can be loaded over the guide wire using the guide wire insertion tool.


Insert the guide wire into the lead by placing the proximal (rigid) end of the guide wire into the distal lead tip using the guide wire insertion tool included in the package (Figure 12).

Note: There may be slight resistance as the guide wire passes through the lead tip seal.

Note: Grasp lead and guide wire insertion tool between your thumb and forefinger to prevent the lead from dislodging from the guide wire insertion tool while you are inserting the guide wire.

Note: Be sure to remove the guide wire insertion tool prior to lead implant.

Figure 12.

- 1 Guide wire
- 2 Lead distal tip

6.10 Placing the lead using a guide wire

Warning: Do not force the lead if significant resistance is encountered during lead passage.

Caution: If the distal end of the guide wire becomes severely kinked or twisted, it may be difficult to withdraw it back through the lead. Therefore, if there is an indication that the distal end of the guide wire has become damaged, or if there is significant resistance in guide wire passage, remove the lead and guide wire. Remove the guide wire from the lead and insert a new guide wire into the lead. Do not use excessive force to retract the guide wire from the lead.

Note: If the lead is not advancing, or if the lead and guide wire seem to be sticking together, there may be thrombus on the guide wire at the lead tip. Remove and inspect the lead and guide wire. Consider using a new guide wire. Reinsert the lead and the guide wire as described in previous steps.

Note: If it is difficult to advance the guide wire around a bend, consider changing the guide wire. More flexible guide wires are recommended for tortuous anatomies. Firmer guide wires are recommended where additional support is needed.

There are many techniques that may be used to advance the lead into a cardiac vein using a stylet. The choice of technique is left to the discretion of the physician. Two possible techniques are described below.

6.10.1 First technique to advance the lead into a cardiac vein using a guide wire

Advance the lead into a cardiac vein:

1. Using venogram images as reference, advance the lead into the coronary sinus. Keep the distal tip of the guide wire beyond the distal tip of the lead during passage and placement (Figure 13).

Figure 13.

2. Rotate the guide wire and advance the guide wire into the cardiac vein. Pass the lead over the guide wire and into the cardiac vein (Figure 14).

Figure 14.

3. Advance the lead further into the cardiac vein by advancing the lead along the guide wire or by advancing the lead and guide wire together.


Note: The proximal electrode should be placed in the selected vein, not in the coronary sinus.

6.10.2 Second technique to advance the lead into a cardiac vein using a guide wire

Advance the lead into a cardiac vein:


1. Using venogram images as reference, advance the lead into the cardiac vein (Figure 15). The curved distal tip of the lead may be used to aid in subselecting cardiac veins within the coronary sinus by slightly retracting the guide wire inside the lead lumen.

Figure 15.

2. If retracted, reinsert the guide wire through the lead tip seal and more distally into the cardiac vein. Advance the lead over the guide wire into the cardiac vein (Figure 16).

Figure 16.

3. Advance the lead further into the cardiac vein by advancing the lead along the guide wire or by advancing the lead and guide wire together.


Note: The proximal electrode should be placed in the selected vein, not in the coronary sinus.

6.11 Fixating the lead

Fixation is achieved the same way regardless of whether a stylet or a guide wire is used.

Using fluoroscopy for guidance, fixate the lead by wedging both curves of the lead tip into the cardiac vein (Figure 17).

Figure 17.


Note: If the selected vein is large, it may be necessary to position the lead in a smaller cardiac vein in order to achieve lead tip fixation. Both electrodes must touch the tissue.

6.12 Taking electrical measurements

The Model 4296 was designed to provide pacing via selectable electrodes. Recommended pacing vectors are shown in Figure 18.

Extended bipolar pacing is available using this lead in combination with a compatible CRT-D system and RV defibrillation lead or with a compatible CRT-P system and RV pacing lead.


Figure 18.

- 1 Extended bipolar, CRT-D — The LV ring is available on specific Medtronic CRT-D devices. Refer to the appropriate CRT-D manual for available LV lead pacing configurations with Medtronic CRT-D devices.
- 2 Unipolar (CRT-P only).
- 3 Extended bipolar, CRT-P — Refer to the appropriate CRT-P manual for available LV lead pacing configurations with Medtronic CRT-P devices.

Warning: The Model 4296 was designed for optimal pacing when used in a unipolar or extended bipolar configuration. The standard bipolar configuration as shown in Figure 19 may result in markedly elevated thresholds or produce anodal stimulation.

Figure 19.

Caution: Prior to taking electrical or defibrillation efficacy measurements, move objects made of conductive materials, such as guide wires or stylets, away from all electrodes.

Note: Initial electrical measurements should be taken with the stylet or guide wire retracted a minimum of 4 cm inside the lead lumen.

Use the following steps to take electrical measurements:

1. Attach a surgical cable to the lead connector pin or ring depending on which electrode you want to use for measurements.
Note: Testing the lead in a unipolar configuration requires the use of an indifferent electrode.
2. Use an implant support instrument, such as a pacing system analyzer, for obtaining electrical measurements. For information on the use of the implant support instrument, consult the product documentation for that device. Satisfactory lead placement is indicated by low stimulation thresholds and adequate sensing of intracardiac signal amplitudes.
 - A low stimulation threshold provides a desirable safety margin, allowing for a possible rise in thresholds that may occur within 2 months following implant.
 - Adequate sensing amplitudes ensure that the lead is properly sensing intrinsic cardiac signals. Minimum signal requirements depend on the sensitivity capabilities of the device. Acceptable acute signal amplitudes for the lead must be greater than the minimum device sensing capabilities. Be sure to include an adequate safety margin to account for lead maturity.

Table 1. Recommended measurements at implant (Assumes 500 Ω resistance)

Measurement recommended	Left ventricle
Maximum acute stimulation thresholds ^a	3.0 V
Minimum acute sensing amplitudes	4.0 mV

^a At a pulse duration setting of 0.5 ms.

3. If electrical measurements do not stabilize to acceptable levels, it may be necessary to reposition the lead and repeat the testing procedure.
4. Check for diaphragmatic stimulation by pacing at 10 V and a pulse width setting greater than 0.5 ms. Then observe for diaphragmatic contracting either by fluoroscopy or direct abdominal palpitation. Further testing may include patient positional changes to simulate upright chronic conditions. If diaphragmatic pacing occurs, reduce the voltage until a diaphragmatic pacing threshold is determined. Diaphragmatic stimulation usually necessitates repositioning of the lead.

6.13 Removing the guide catheter from the lead

Note: If an Attain Hybrid guide wire was used for placement of the lead, go to Step 2. It is not necessary to replace the Attain Hybrid guide wire with a stylet for catheter removal.

Once the lead is in the final position, remove the guide catheter from the lead:

1. If used, remove the guide wire and guide wire insertion tool. Replace the guide wire with a straight stylet (downsized knob). Insert the straight stylet into the lead to the mid-coronary sinus.
2. Remove the guide catheter from the lead. See the delivery system product documentation for details.
Note: For Medtronic slittable delivery systems, use a slitter compatible with a 1.75 mm (5.3 French) lead body.
3. Carefully and completely remove the stylet or Attain Hybrid guide wire. When removing the stylet, grip the lead firmly just distal of the connector pin; this will help prevent possible lead tip dislodgement.
4. Repeat the electrical measurements. See Section 6.12, "Taking electrical measurements", page 12.

6.14 Anchoring the lead

Caution: Use care when anchoring the lead.

- Use only nonabsorbable sutures to anchor the lead.
- Do not attempt to remove or cut the anchoring sleeve.
- Do not use the anchoring sleeve tabs for suturing.
- During lead anchoring, take care to avoid dislodging the lead tip.
- Do not secure sutures so tightly that they damage the vein, lead, or anchoring sleeve (Figure 20).
- Do not tie a suture directly to the lead body (Figure 20).


Figure 20.

Anchor the lead using all 3 grooves:

1. Position the anchoring sleeve against or near the vein.
2. Secure the anchoring sleeve to the lead body by tying a suture firmly in each of the 3 grooves (Figure 21).

Figure 21.

3. Use at least 1 additional suture in 1 of the grooves to secure the anchoring sleeve and lead body to the fascia.

6.15 Connecting the lead

Caution: Always remove the stylet and stylet guide before connecting the lead to the device. Failure to remove the stylet may result in lead failure.

Connect the lead to an implantable device.


1. Insert the lead connector into the connector block. Consult the product documentation packaged with the implantable device for instructions on proper lead connections.
2. Take electrical measurements through the device.

6.16 Placing the device and lead into the pocket

Caution: Use care when placing the device and leads into the pocket.

- Ensure that the leads do not leave the device at an acute angle.
- Do not grip the lead or device with surgical instruments.
- Do not coil the lead. Coiling the lead can twist the lead body and may result in lead dislodgement (Figure 22).

Figure 22.

Use the following steps to place the device and leads into the pocket:

English 13

1. To prevent undesirable twisting of the lead body, rotate the device to loosely wrap the excess lead length (Figure 23).

Figure 23.

2. Insert the device and leads into the pocket.
3. Before closing the pocket, verify sensing, pacing, cardioversion, and defibrillation efficacy.

6.17 Post-implant evaluation

After implant, monitor the patient's electrocardiogram until the patient is discharged. If a lead dislodges, it usually occurs during the immediate postoperative period.

Recommendations for verifying that the lead is properly positioned include x-rays and pacing thresholds taken at pre-hospital discharge, 3 months after implant, and every 6 months thereafter.

In the event of a patient death, explant all implanted leads and devices and return them to Medtronic with a completed Product Information Report form.

7 Specifications

7.1 Specifications (nominal)


Parameter	4296
Type	Dual electrode
Chamber paced	Left ventricle
Length	78 and 88 cm
Connector	IS-1 BI
Material	Conductor: 25% Ag-core-MP35N
	Insulators: Polyurethane (outer) SI-polyimide (SI-PI) (inner) ^a
	Ring electrode: Platinum iridium alloy with titanium nitride coating
	Tip electrode: Platinum iridium alloy with titanium nitride coating
	Connector pin: Stainless steel
	Connector ring: Stainless steel
	Molded tip seal: Silicone rubber (with steroid)

Parameter	4296
Tip electrode configuration	Tapered, annular, titanium-nitride-coated, steroid-eluting
Ring electrode configuration	Radiused, titanium-nitride-coated, steroid-eluting
Diameter	Lead body: 1.75 mm (5.3 French)
	Tip electrode: 1.53 mm (4.6 French)
	Ring electrode: 1.70 mm (5.1 French)
	1.90 mm (5.7 French)
Medtronic delivery system (recommended inner diameter)	
Diagnostic guide wire (recommended diameter)	0.36 mm to 0.46 mm
Electrode surface area	Tip: 5.8 mm ²
	Ring: 5.8 mm ²
Conductor resistance	Tip: 38 ± 10 Ω (78 cm)
	43 ± 10 Ω (88 cm)
	Ring: 35 ± 10 Ω (78 cm)
	40 ± 10 Ω (88 cm)
Distance between electrodes	21 mm
Steroid	Dexamethasone acetate
Combined total of steroid	<1.0 mg

^a Technology developed by NASA.

7.2 Specifications drawing (nominal)

Figure 24.

- 1 Proximal electrode: nominally 5.8 mm² geometric pacing surface area
- 2 Tip electrode: nominally 5.8 mm² geometric pacing surface area
- 3 Anchoring sleeve
- 4 IS-1 BI connector

8 CE mark of conformity

The following CE mark of conformity applies to the Model 4296 lead.

CE0123
2009

9 Explanation of symbols on package labeling

Refer to the package labels to see which symbols apply to this product.

Table 2. Explanation of symbols on package labeling

Symbol	Explanation
CE0123	Conformité Européenne (European Conformity). This symbol means that the device fully complies with applicable European Union acts.
	Do not use if package is damaged
	Do not reuse
	Upper limit of temperature
STERILE EO	Sterilized using ethylene oxide
	Consult instructions for use
	Consult instructions for use at this website
	Date of manufacture
	Manufacturer
	Use by
REF	Reorder number
EC REP	Authorized representative in the European Community
SN	Serial number

Table 2. Explanation of symbols on package labeling (continued)

Symbol	Explanation
	Open here
	Package contents
	Product documentation
	Caution
	Accessories
	Inner diameter
	Lead length
	Catheter delivered OTW
	Catheter
	Lead
	Transvenous cardiac vein pacing lead
	Transvenous ventricular lead
	Transvenous atrial lead
	Transvenous lead with one defibrillation electrode
	Transvenous lead with two defibrillation electrodes
	Pace
	Sense

Table 2. Explanation of symbols on package labeling (continued)

Symbol	Explanation
	Defibrillation
	Extendable and retractable screw-in
	Tined
	Steroid eluting
	Lead introducer
	Lead introducer with guide wire
	J-shaped
	SureScan logo
	MR Conditional

10 Medtronic disclaimer of warranty

For complete disclaimer of warranty information, see the accompanying disclaimer of warranty document.

11 Service

Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed on the back cover.

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
Declaration Letter

Medtronic Inc., a corporation organized under the laws of the United States, having its principal place of business at 710 Medtronic Parkway, Minneapolis, MN 55432, USA, hereby or affirm that the attached documentation, IFU for MRI model 4396, is true and unaltered copies of the original document.


Name: Shellee Lamb

Title: Sr. Regulatory Affairs Specialist

Oct 4, 2019

Date

Cindi A Digelbach
10/4/19

Medtronic

ATTAIN ABILITY™ STRAIGHT MRI SURESCAN™ 4396

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Contents

1	Description	3
2	Indications	4
3	Contraindications	4
4	Warnings and precautions	4
5	Potential adverse events	6
6	Implant procedure	7
7	Specifications	13
8	CE mark of conformity	14
9	Explanation of symbols on package labeling	14
10	Medtronic disclaimer of warranty	15
11	Service	15

1 Description

The Medtronic Attain Ability Straight MRI SureScan 4396 steroid eluting, dual electrode, transvenous, over the wire, cardiac vein pacing lead is designed for pacing and sensing via a cardiac vein. The lead has been tested for use in the Magnetic Resonance Imaging (MRI) environment. All lead lengths for this lead model are MR Conditional. This lead contains 2 electrodes: the distal electrode, positioned at the distal tip, and the proximal electrode, positioned 21 mm from the tip. See Section 7.2, "Specifications drawing (nominal)", page 14 for a lead drawing.

The distal tip electrode allows a guide wire to pass through to aid in cardiac vein selection. The distal tip electrode also contains a dual-function steroid-loaded silicone rubber membrane MCRD¹. This membrane seals the lead inner lumen to reduce blood ingress, and the MCRD elutes steroid to reduce inflammatory response within the cardiac vein. The proximal electrode is proximal to the tines. It also contains an MCRD for elution of steroid. The outer insulation of the lead is polyurethane and the inner insulation is SI-polyimide (SI-PI)². The SI-PI is applied as a coating to the conductor wire prior to coiling. The IS-1 connector pin has a drilled inner lumen to improve guide wire passage³.

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The Model 4396 lead can be positioned with the aid of a guide wire or with a stylet. If a stylet is used, use only the stylets packaged with the lead or in a stylet kit (downsized knob). Always use a stylet that is 3 cm shorter than the lead length listed on the IS-1 connector label.

Note: To implant the Model 4396 lead in a cardiac vein, a compatible delivery system is required. A compatible delivery system includes a guide catheter and either a hemostasis valve or an introducer valve that can be removed or that allows passage over the IS-1 connector. Contact a Medtronic representative for further information regarding compatible delivery systems.

1.1 Medtronic SureScan system

The Model 4396 lead is part of the Medtronic SureScan system. The SureScan system includes a Medtronic SureScan device connected to Medtronic SureScan leads. Labeling for SureScan system components displays the SureScan logo and the MR Conditional symbol.

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The Model 4396 lead is contraindicated as follows:

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Handling steroid MCRDs⁴ – Avoid reducing the amount of steroid available before implanting the lead. Reducing the available amount of steroid may adversely affect low-threshold performance.

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Bipolar pacing – This lead has not been designed for pacing in the bipolar configuration. It utilizes 2 equal-sized electrodes. Use

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of this lead in a bipolar configuration may result in higher pacing thresholds or anodal stimulation.

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- Rust stylets are not recommended with this lead due to the risk of conductor coil or insulation perforation.
- If a stylet is used for lead positioning, use only the stylets packaged with the lead or in a stylet kit (downsized knob). Other stylets may extend beyond the lead tip causing lead tip seal damage or patient injury.
- If the lead is damaged, do not implant it. Return the lead to a Medtronic representative.
- Protect the lead from materials that shed particles such as lint and dust. Lead insulators attract these particles.
- Handle the lead with sterile surgical gloves that have been rinsed in sterile water or a comparable substance.
- Do not severely bend, kink, or stretch the lead.
- Do not use surgical instruments to grasp the lead or connector pin.

- Do not immerse leads in mineral oil, silicone oil, or any other liquid, except blood, at the time of implant.
- Use an anchoring sleeve with all leads. Ensure that the anchoring sleeve is positioned close to the lead connector pin, to prevent inadvertent passage of the sleeve into the vein. If wiping the lead is necessary prior to insertion, ensure that the anchoring sleeve remains in position.
- Do not force the guide catheter or leads if significant resistance is encountered. Use of guide catheters or leads may cause trauma to the heart.

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Diathermy treatment (including therapeutic ultrasound) – Diathermy is a treatment that involves the therapeutic heating of body tissues. Diathermy treatments include high frequency, short wave, microwave, and therapeutic ultrasound. Except for therapeutic ultrasound, do not use diathermy treatments on cardiac device patients. Diathermy treatments may result in serious injury or damage to an implanted device and lead system. Therapeutic ultrasound (including physiotherapy, high intensity therapeutic ultrasound, and high intensity focused ultrasound), is the use of ultrasound at higher energies than diagnostic ultrasound to bring heat or agitation into the body. Therapeutic ultrasound is acceptable if treatment is performed with a minimum separation distance of 15 cm between the applicator and the

English 5

implanted device and lead system, as long as the ultrasonic beam is pointing away from the device and lead system.

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- air embolism
- cardiac dissection
- cardiac perforation
- cardiac tamponade
- coronary sinus dissection
- death
- endocarditis and pericarditis
- erosion through the skin
- extracardiac muscle or nerve stimulation
- fibrillation or other arrhythmias
- heart block
- heart wall or vein wall rupture
- hematoma/seroma
- infection
- myocardial irritability
- myopotential sensing
- pericardial effusion
- pericardial rub
- pneumothorax
- body rejection phenomena including local tissue reaction
- threshold elevation
- thrombosis
- thrombotic embolism
- valve damage (particularly in fragile hearts)

Potential adverse events related to the lead and the programmed parameters include, but are not limited to, the following events:

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Implant techniques that may damage the lead	Possible effects on the lead	Corrective action to be considered
Forcing the lead through the introducer/delivery system	Electrode, conductor coil, and/or insulation damage	Replace the lead.
Use of too medial of an approach with venous introducer resulting in clavicle and first rib binding	Conductor coil fracture, insulation damage	Replace the lead.
Using too stiff a stylet	Conductor coil/insulation perforation	Replace the lead.
Puncturing the pericardium and/or tendon when using subclavian introducer approach	Conductor coil fracture, insulation damage	Replace the lead.
Advancing the lead through the non-coronary central access veins without the stylet or guide wire fully inserted	Tip distortion and/or insulation perforation	Replace the lead.
Inserting the proximal end of the guide wire through the lead tip seal without using the guide wire insertion tool	Lead tip seal damage and/or conductor coil/insulation damage	Replace the lead.
Advancing a stylet tip beyond the distal end of the lead tip seal	Lead tip seal damage	Replace the lead.

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Note: To implant the Model 4396 lead in a cardiac vein, a compatible delivery system is required, such as a Medtronic delivery system. A compatible delivery system includes a guide catheter and a hemostasis valve or introducer valve which allows passage through or removal from an IS-1 connector. Contact a Medtronic representative for further information regarding compatible delivery systems.

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1. Insert a straight stylet into the lead.
Note: When the stylet is fully inserted, the distal tip of the stylet does not reach the distal tip of the lead.
2. Insert the lead into the delivery system. See the delivery system product documentation for additional information.

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
The Model 4396 lead can be placed with the aid of a stylet or with a guide wire.

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Stylet delivery – If the patient's anatomy features gentle vein angulation off of the coronary sinus and the cardiac vein branch is not tortuous (Figure 2), use a stylet for lead delivery.


Figure 2.

- 1 Coronary sinus
- 2 Cardiac vein

Guide wire delivery – If the patient's anatomy features acute vein angulation off of the coronary sinus and the cardiac vein branch is tortuous (Figure 3), use a guide wire for lead delivery.

Figure 3.

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6.8 Placing the lead using a stylet

A stylet may be used to place the lead into a cardiac vein. The physician must shape the stylet to place the lead into an appropriate small, tortuous cardiac vein.

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
Caution: To avoid damage to the stylet, do not use a sharp object to impart a curve to the distal end of the stylet.

Note: If a slight curve is needed for the stylet, use only a smooth object to impart a curve to the distal portion of a stylet (Figure 4).

Figure 4.

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
There are many techniques that may be used to advance the lead into a cardiac vein using a guide wire. The choice of technique is left to the discretion of the physician. One possible technique is described below.

6.8.1 Using a stylet to advance the lead into the selected cardiac vein


1. Advance the lead up to, but not past, the cardiac vein ostium (Figure 5).

Figure 5.

2. Use the following actions to direct the lead tip into the cardiac vein:
 - Remove the straight stylet.
 - Manually shape the stylet into a curve.
 - Insert the curved stylet into the lead.
 - Rotate the lead body and stylet together.
3. Advance the lead into the selected cardiac vein (Figure 6).

Figure 6.

4. Advance the lead further into the cardiac vein by advancing the lead along the stylet or by advancing the lead and stylet together.

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6.9 Placing the lead using a guide wire

Warning: Do not insert the proximal end of the guide wire through the lead tip seal without using the guide wire insertion tool.

Inserting the guide wire without the guide wire insertion tool could cause damage to the lead tip seal or to the conductor core or insulation.

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1. Select a guide wire. A more flexible guide wire is recommended if the patient has tortuous anatomy. If additional support is needed, use a firmer guide wire.

Note: Guide wire firmness is determined by 2 factors: guide wire diameter and guide wire core design. A larger guide wire diameter may be firmer than a smaller guide wire diameter. However, guide wires with the same diameter may have different degrees of firmness. A "J" shaped guide wire or angled tip guide wire are recommended to aid in sub-selecting the cardiac vein and where more steerability is desired.

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
2. Insert the guide wire into the lead by placing the distal (flexible) end of the guide wire into the lead connector pin using the guide wire insertion tool included in the package (Figure 7).

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
5. Position the guide wire steering handle (Figure 8) using the following steps :
 - a. Advance the guide wire steering handle over the proximal (rigid) end of the guide wire.
 - b. Tighten the guide wire steering handle onto the guide wire near the lead connector pin.

Figure 8.

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
Insert the guide wire into the lead by placing the proximal (rigid) end of the guide wire into the distal lead tip using the guide wire insertion tool included in the package (Figure 9).

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Figure 9.

- 1 Guide wire
- 2 Lead distal tip

6.9.2 Using a guide wire to advance the lead into the selected cardiac vein

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Caution: If the distal end of the guide wire becomes severely kinked or twisted, it may be difficult to withdraw it back through the lead. Therefore, if there is an indication that the distal end of the guide wire has become damaged, or if there is significant resistance in guide wire passage, remove the lead and guide wire. Remove the guide wire from the lead and insert a new guide wire.

into the lead. Do not use excessive force to retract the guide wire from the lead.

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There are many techniques that may be used to advance the lead into a cardiac vein using a guide wire. The choice of techniques is left to the discretion of the physician. One possible technique is described below.


1. Using venogram images as reference, advance the lead into the coronary sinus. Keep the distal tip of the guide wire beyond the distal tip of the lead during passage and placement (Figure 10).

Figure 10.

2. Rotate the guide wire and advance the guide wire into the cardiac vein. Pass the lead over the guide wire and into the cardiac vein (Figure 11).

Figure 11.

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
Using fluoroscopy for guidance, fixate the lead into the cardiac vein (Figure 12).

Note: If the selected vein is large, it may be necessary to position the lead in a smaller, more tortuous cardiac vein in order to achieve lead tip fixation. Both electrodes must touch the tissue.

Figure 12.


6.11 Taking electrical measurements

The Model 4396 was designed to provide pacing via selectable electrodes. Recommended pacing vectors are shown in Figure 13.

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
Warning: The Model 4396 was designed for optimal pacing when used in a unipolar or extended bipolar configuration. The standard bipolar configuration as shown in Figure 14 may result in markedly elevated thresholds or produce anodal stimulation.

Figure 14.

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^a At a pulse duration setting of 0.5 ms.

- If electrical measurements do not stabilize to acceptable levels, it may be necessary to reposition the lead and repeat the testing procedure.
- Check for diaphragmatic stimulation by pacing at 10 V and a pulse width setting greater than 0.5 ms. Then observe for diaphragmatic contracting either by fluoroscopy or direct abdominal palpitation. Further testing may include patient positional changes to simulate upright chronic conditions. If diaphragmatic pacing occurs, reduce the voltage until a diaphragmatic pacing threshold is determined. Diaphragmatic stimulation usually necessitates repositioning of the lead.

6.12 Removing the guide catheter from the lead

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- If used, remove the guide wire and guide wire insertion tool. Replace the guide wire with a straight stylet (downsized knob). Insert the straight stylet into the lead to the mid-coronary sinus.
- Remove the guide catheter from the lead. See the delivery system product documentation for details.
Note: For Medtronic slittable delivery systems, use a slitter compatible with a 1.3 mm (4 French) lead body.
- Carefully and completely remove the stylet or Attain Hybrid guide wire. When removing the stylet, grip the lead firmly just distal of the connector pin; this will help prevent possible lead tip dislodgement.
- Repeat the electrical measurements. See Section 6.11, "Taking electrical measurements", page 11.

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
- Use only nonabsorbable sutures to anchor the lead.
- Do not attempt to remove or cut the anchoring sleeve.
- Do not use the anchoring sleeve tabs for suturing.
- During lead anchoring, take care to avoid dislodging the lead tip.
- Do not secure sutures so tightly that they damage the vein, lead, or anchoring sleeve (Figure 15).
- Do not tie a suture directly to the lead body (Figure 15).

Figure 15.

Anchor the lead using all 3 grooves:

- Position the anchoring sleeve against or near the vein.
- Secure the anchoring sleeve to the lead body by tying a suture firmly in each of the 3 grooves (Figure 16).

Figure 16.

- Use at least 1 additional suture in 1 of the grooves to secure the anchoring sleeve and lead body to the fascia.

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
- Insert the lead connector into the connector block. Consult the product documentation packaged with the implantable device for instructions on proper lead connections.
- Take electrical measurements through the device.

6.15 Placing the device and lead into the pocket

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- Ensure that the leads do not leave the device at an acute angle.
- Do not grip the lead or device with surgical instruments.
- Do not coil the lead. Coiling the lead can twist the lead body and may result in lead dislodgement (Figure 17).

Figure 17.

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1. To prevent undesirable twisting of the lead body, rotate the device to loosely wrap the excess lead length (Figure 18).

Figure 18.

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7 Specifications

7.1 Specifications (nominal)

Parameter	Model 4396
Type	Dual electrode
Chamber paced	Left ventricle
Length	78 and 88 cm
Connector	IS-1 BI
Material	Conductor: 25% Ag-core-MP35N
	Insulators: Polyurethane (outer) SI-polyimide (SI-PI) (inner) ^a
	Ring electrode: Platinum iridium alloy with titanium nitride coating
	Tip electrode: Platinum iridium alloy with titanium nitride coating
	Connector pin: Stainless steel
	Connector ring: Stainless steel
	Molded tip seal: Silicone rubber (with steroid)
Tip electrode configuration	Tapered, annular, titanium-nitride-coated, steroid-eluting


Parameter	Model 4396
Ring electrode configuration	Radiused, titanium-nitride-coated, steroid-eluting
Diameter	Lead body: 1.3 mm (4 French) Tip electrode: 1.53 mm (4.6 French) Ring electrode: 1.70 mm (5.1 French)
Medtronic delivery system (recommended inner diameter)	1.90 mm (5.7 French)
Diagnostic guide wire (recommended diameter)	0.36 mm to 0.46 mm
Electrode surface area	Tip: 5.8 mm ² Ring: 5.8 mm ²
Conductor resistance	Tip: 38 ± 10 Ω (78 cm) 43 ± 10 Ω (88 cm) Ring: 35 ± 10 Ω (78 cm) 40 ± 10 Ω (88 cm)
Distance between electrodes	21 mm
Steroid	Dexamethasone acetate
Combined total of steroid	<1.0 mg

^a Technology developed by NASA.

English 13

7.2 Specifications drawing (nominal)

Figure 19.

- 1 Tip electrode: nominally 5.8 mm² geometric pacing surface area
- 2 Proximal ring electrode: nominally 5.8 mm² geometric pacing surface area
- 3 Anchoring sleeve
- 4 IS-1 BI connector

8 CE mark of conformity

The following CE mark of conformity applies to the Model 4396 lead.

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Symbol	Explanation
CE0123	Conformité Européenne (European Conformity) This symbol means that the device fully complies with applicable European Union acts.
	Do not use if package is damaged
	Do not reuse
	Upper limit of temperature
STERILE EO	Sterilized using ethylene oxide
	Consult instructions for use
	Consult instructions for use at this website
	Date of manufacture
	Manufacturer
	Use by
REF	Reorder number
EC REP	Authorized representative in the European Community
SN	Serial number

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© 2015 Medtronic, Inc.