**Legal Manufacturer:**

Siemens Healthcare GmbH  
Henkestr. 127, 91052 Erlangen  
Germany

***User Manual Supplement for Russia***

*The MR-PET system Biograph mMR with accessories*

*M. Velasco*

---

Mayra Alejandra Velasco Alvarado  
Regulatory Affairs  
Siemens Healthcare GmbH, Magnetic Resonance (MR)

## **1 The name of the Medical Device**

The MR-PET system Biograph mMR with accessories.

## **2 Manufacturing information of the Medical Device**

### **2.1 Legal manufacturer**

Name: Siemens Healthcare GmbH  
Address: Henkestr. 127, 91052 Erlangen, Germany  
Phone: +49 9131 84-0  
Web: [www.siemens-healthineers.com](http://www.siemens-healthineers.com)

### **2.2 Manufacturing site**

Name: Siemens Healthcare GmbH, Magnetic Resonance (MR)  
Address: Henkestr. 127, 91052 Erlangen, Germany  
Phone: +49 9131 84-0  
Web: [www.siemens-healthineers.com](http://www.siemens-healthineers.com)

### **2.3 Local representative in Russia**

Name: Siemens Healthcare LLC  
Address: Dubininskaja street 96, 115093 Moscow, Russia  
Tel.: +7 495 737 12 52  
Web: [www.siemens-healthineers.com/ru/](http://www.siemens-healthineers.com/ru/)

## **3 Designation of the medical device**

### **3.1 Indications and intended use**

The Siemens MR-PET system combines magnetic resonance diagnostic devices (MRDD) and Positron Emission Tomography (PET) scanners that provide registration and fusion of high resolution physiologic and anatomic information, acquired simultaneously and isocentrically. The combined system maintains independent functionality of the MR and PET devices, allowing for single modality MR and / or PET imaging.

These systems are intended to be utilized by appropriately trained health care professionals to aid in the detection, localization, and diagnosis of diseases and disorders.

The MR is intended to produce transverse, sagittal, coronal and oblique cross-sectional MR-images, spectroscopic images and/or spectra, and displays the internal structure and/or function of the human body

Other physical parameters derived from the images and/or spectra may also be produced. Depending on the region of interest, approved contrast agents may be used, as described in their labeling. This system may also be used for imaging during interventional procedures when performed with MR compatible devices, such as MR Safe biopsy needles.

The PET images and measures the distribution of PET radiopharmaceuticals in humans to aid the physician in determining various metabolic (molecular) and physiologic functions within the human body for evaluation of diseases and disorders such as, but not limited to, cardiovascular disease, neurological disorders and cancer. The combined system utilizes the MR for radiation-free attenuation correction maps for PET studies. The system provides inherent anatomical reference for the fused PET and MR images due to precisely aligned MR and PET image coordinate systems.

The MR-PET system is not a device with measuring function as defined in the Medical Device Directive (MDD). Quantitative measured values obtained are for informational purposes and cannot be used as the only basis for diagnosis.

For the USA only: Federal law restricts this device to sale, distribution and use by or on the order of a physician.

Your system is a medical device for human use only!

### 3.1.1 MR mammography

MR mammography is the most sensitive method for detecting and staging invasive breast carcinomas. Benign lesions can be differentiated from malignant lesions using a combination of anatomical, functional, and metabolic information.

A special application is used in areas around prosthetic implants.

#### Indications:

- Verification of invasive carcinoma.
- Prosthetic implants (complications, exclusion of tumors).
- Evaluation of benignity and malignancy in cases of structural thickening.
- Preoperative determination of multicentricity or multifocality for breast-preserving therapy.
- Distinction between mastopathies and malignant tumor formations.
- Differentiation of scar tissue.
- Follow-up control after tumor surgery.

### 3.2 Contraindications

This chapter contains very important information about MR safety: In principle, a qualified physician must evaluate the risk/benefit ratio of the MR examination for every patient. To date, there is no scientific proof that MR examinations are harmless for pregnant women, the unborn (embryos or fetuses) and children under two years of age.

For MR-PET examinations, also observe the specific PET-related safety advices.

In general, an MR examination is contraindicated for patients with electronic or electronically conductive implants or metals, especially those containing ferromagnetic foreign matter.

Typical contraindications for MR examinations are:

- Electronic implants: e.g. pacemakers, stimulators, insulin pumps
- Artificial heart valves, aneurysm clips
- Metal splinters in the eye (danger of retinal detachment)
- Artificial anus (anus praeter) with magnetic closure
- Transdermal drug patches with metallic backings
- Electrically conductive implants and prostheses
- Metallic spirals for contraception (IUDs = Intrauterine Devices)
- Transdermal or other similar implants (for example, body piercings as well as magnetic piercings)

#### **WARNING**

Electronic and/or electrically conductive implants and magnetizable inclusions in static and low-frequency magnetic fields and RF fields!

##### **Risk of patient death/patient injury**

- ◆ Ask the patient about implants and inclusions.
- ◆ Do not perform MR examinations on patients with electronic or electrically conductive implants and magnetizable inclusions.
- ◆ Ensure that patients wearing such implants and/or inclusions remain outside the exclusion zone (0.5 m line).

Exceptions: Certain implantable medical devices have been cleared, approved and/or licensed by the Competent Governmental Authorities and/or labeled by the device manufacturer as "MR Conditional". For such implantable medical devices, the previously mentioned list of general contraindications and the warning may not be applicable in its entirety.

It is the responsibility of the device manufacturer to declare an implantable medical device as MR Conditional, appropriate and to define the conditions (constraints) for safe MR scanning. The MR operator must be aware of any such conditions for MR scanning. It is the obligation of the MR operator to assure that these conditions are strictly adhered to. To obtain these specific conditions the MR operator may refer to the labeling of the implantable medical device or contact the device manufacturer. Siemens MR does not assume responsibility or liability for the operation of the MR system with any implantable medical device. Especially, Siemens MR is not responsible for controlling technical parameters of the MR system other than those defined by the normal operating mode, the first level controlled operating mode and the data provided in the system owner manual, such as spatial gradient field plots.

**CAUTION**

Eddy currents induced by low-frequency magnetic fields!

**Patient burns**

- ◆ Do not examine patients with electrically conducting implants or prostheses.

**CAUTION**

Electrically conducting objects!

**Injury to patient due to warming**

**Incorrect diagnosis due to artifacts**

- ◆ Request that the patient removes all electrically conducting objects, e.g. necklaces, rings, braces, rubber bands for long hair, piercings as well as jewelry.
- ◆ Request that the patient removes all clothing including electrically conducting material, for example, bras, metallic appliqués or woven metallic yarns.
- ◆ Inform patients that eyeliners and tattoos may contain ingredients causing artifacts or skin irritations during MR examinations. In some cases, patients have been burned.
- ◆ To prevent injuries, instruct patients to remove makeup prior to the examination.
- ◆ Instruct patients to seek medical attention in case of discomfort during or following the MR examination.

### 3.2.1 MR mammography

Contra-indications:

- Evaluation of microcalcification.
- Inflammatory changes.
- Mammary secretions.
- Known mastopathies.

### 3.3 Possible side effects

Possible undesirable side effects for MR are dizziness, heating, claustrophobia and nerve stimulation.

### 3.4 Requirements to be met by MR mammography

- Full coverage of the breast.
- High accuracy for differentiating breast tissue, especially for patients with dense breasts or breast implants.

### 3.5 Typical workflow of routine clinical breast MR examinations:

- Characterization of tumor cellularity with diffusion-weighted imaging.
- Imaging of kinetic patterns (dynamic measurements).
- High-resolution 3D MR imaging of morphology.
- Characterization of tumor metabolism with 1H MR spectroscopy.
- With prosthetic implants: Display of silicone images.

### 3.6 Physiological effects

Due to the presence of alternating electromagnetic fields, patients may experience various physiological effects during MR measurements:

- Peripheral nerve stimulation through low-frequency fields of the gradient coils
- Warming of body tissue through RF fields of the RF transmitter coil

These physiological effects can be evaluated by the technical quantities dB/dt (stimulations) and SAR (warming) respectively.

It is generally accepted that no published evidence supporting the occurrence of cumulative and/or long-term effects after exposure to EMF emitted by the MR equipment exists.

### 3.7 Information about PET safety

This chapter includes the description of specific PET hazards and risks which occur during an MR-PET operation.

In addition to the PET-specific advices, always observe the MR related safety advice notices.

A PET system does not present permanently existing hazards and risks. However, the user must pay particular attention to the following topics:

- Radioactive material/tracer
- Attenuation correction and imaging artifacts

### 3.8 Authorized operating personnel

The MR-PET system must be operated according to the intended use and only by qualified persons with the necessary knowledge in accordance with country-specific regulations, e.g. physicians, trained radiological technicians or technologists, subsequent to the necessary user training.

This user training must include basics in MR technology, PET technology as well as safe handling of MR systems. The user must be familiar with potential hazard and safety guidelines the same way the user is familiar with emergency and rescue scenarios. In addition, the user has to have read and understood the contents of the operator manual.

Please contact Siemens Service for more information on available training options and suggested duration and frequency of such training.

### 3.9 Mode of application (Medical Device application data)

1. Preparing an examination
  - Patient registration
  - Positioning the patient on the patient table
  - Positioning the local coils and cushions
  - Table positioning
2. Performing an examination
  - Measuring of the Localizer and parallel loading of previous reference images
  - Planning the examination volume
  - Scanning
  - Viewing and editing images during examinations
3. Basic MR viewing, post-processing, and filming
  - Viewing

- Mean Curve
- Argus Viewer
- Dynamic Analysis
- 3D
- Filming

### 3.10 Ambient conditions

|                       | Control room               | Electronics room           | Examination room           |
|-----------------------|----------------------------|----------------------------|----------------------------|
| Temperature           | 15-30 °C                   | 15-30 °C                   | 18-22 °C                   |
| Relative air humidity | 40 to 60%                  | 40 to 80%                  | 40 to 60%                  |
| Absolute air humidity | <11.0 g/kg                 | <11.0 g/kg                 | <11.0 g/kg                 |
| Air pressure          | 700-1060 hPa <sup>1)</sup> | 700-1060 hPa <sup>1)</sup> | 700-1060 hPa <sup>1)</sup> |

\* - 700 hPa is equivalent to a height of about 3000m above sea level.

## 4 Classification of the Medical Device

### 4.1 Type of the medical device in accordance with the nomenclature for classification

#### 4.1.1 Europe (CE)

|                                                            |                                                                                                                                                                                                                          |
|------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Kind of product:                                           | Active medical device for diagnosis                                                                                                                                                                                      |
| Classification rules according to MDD 93/42/EEC, Annex IX: | Rule 10                                                                                                                                                                                                                  |
| Classification according to Annex IX:                      | Ila                                                                                                                                                                                                                      |
| Classification criteria:                                   | Active devices intended for diagnosis are in Class Ila if they are intended to supply energy which will be absorbed by the human body, except for devices used to illuminate the patient's body, in the visible spectrum |

#### 4.1.2 Russia

|                                                                       |        |
|-----------------------------------------------------------------------|--------|
| Classification in accordance with Russia Nomenclature classification: | 221690 |
| Classification in accordance with potential risk of using:            | 2a     |

### 4.2 Software classification according to IEC 62304:2006

The Medical Device's Software Safety Class according to IEC 62304 is Class B.

### 4.3 Software version

Brand name: syngo MR

Software full version: E11P.

Software release version: E11.

Date of latest release: 02/2020.

#### 4.4 Class of the Medical Device in accordance with IEC 60601-1

According to IEC 60601-1, the Medical Device is the **Class 1**.

#### 4.5 Type on contact with the human body

The Medical Device is non-sterile, nonexpendable.

#### 4.6 Laser light localizer

Laser light Class **2M** according to IEC 60825.

#### 4.7 IP protection class

IP protection class according to IEC 60529:

**IPX0** for the system.

#### 4.8 Protection class

Protection class I according to IEC 60601-1: system.

Protection class II according to IEC 60601-1: mMR Breast Coil; Tx/Rx CP Extremity Coil; Shoulder Array Coil (small and large).

#### 4.9 Applied parts

According to IEC 60601-1 Cl. 7.2.10 and 8.3 the following parts are applied parts of the following types:

|                                                                                                                                                                                                                                                                                                                                                                          |                      |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------|
| Patient table MR Tim.                                                                                                                                                                                                                                                                                                                                                    | Type B applied part  |
| Local coils:<br>mMR Head/Neck, A Tim coil (basic and additional)<br>mMR Head 32, A 3T Tim Coil<br>mMR Spine, A Tim coil (basic and additional)<br>mMR Body, A Tim coil (basic and additional)<br>3T Flex large 4 (basic and additional)<br>3T Flex small 4 (basic and additional)<br>mMR Breast Coil<br>Tx/Rx CP Extremity Coil<br>Shoulder Array Coil (small and large) | Type BF applied part |
| Body coil                                                                                                                                                                                                                                                                                                                                                                | Type B applied part  |
| ECG/Pulse module                                                                                                                                                                                                                                                                                                                                                         | Type BF applied part |

#### 4.10 Explosion protection

The MR system is not intended for operation in areas prone to explosion (e.g., highly flammable mixtures of anesthesia gases with air or oxygen or nitrous oxide)

#### 4.11 Disinfection receptors/system components

Disinfectants without alcohol, ether.

#### 4.12 Degree of safety in the presence of a flammable anesthetic's mixture with air or with oxygen or with nitrous oxide

No AP or APG category equipment

#### 4.13 Operating mode

Continuous operation with short-term load


## 5 Technical Data: Specifications and Characteristics

### 5.1 Photos and line illustrations of the device

Array Coil

; types:

## 5.2 The list of components and accessories

| Name of the component and its description |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Basic block</b>                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| 1.                                        | Biograph mMR – System<br>Whole-body superconductive 3T magnet with active shielding (AS) technology with counter coils.                                                                                                                                                                                                                                                                                                                                                                                                 |
| 2.                                        | PET detectors:<br>material: LSO<br>Crystal element dimension: 4 × 4 × 20 mm<br>Crystal elements per block: 64<br>Avalanche Photodiodes (APDs): 9 per block<br>Detector ring diameter: 656 mm<br>Blocks per detector ring: 56<br>Transaxial FOV: 594 mm<br>Plane spacing: 2.0 mm<br>Crystal elements per ring: 448<br>Qty detector element rings: 64<br>Total detector blocks: 448<br>Total # of crystal elements: 28 672<br>Image planes: 127                                                                           |
| 3.                                        | Tim [102x18] MQ Gradients<br>Tim [102x32] MQ Gradients<br>MQ Gradient System with noise reduction features.<br>Integrated measures comprehensively reduce acoustic noise without compromising strongest gradient performance.                                                                                                                                                                                                                                                                                           |
| 4.                                        | Tim [102x18] MQ Gradients<br>Tim [102x32] MQ Gradients<br>Total imaging matrix with 102 channels (coil elements) that can be connected simultaneously, as well as 18 (32) independent RF channels that can be used simultaneously in one single scan and in one single FOV, each generating an independent partial image. It is for demanding high-end applications and optimized throughput. Tim [102x18] has flexibility in Parallel Imaging. Maximum SNR is ensured through the matrix coil and mMR coil technology. |
| 5.                                        | Pure White Design #mMR<br>The color of the main face plate cover with integrated control panel and table display is brilliant white surrounded by a brilliant silver trim. The asymmetrical deco area on the left side is colored white matte and also with a brilliant surrounding silver trim.                                                                                                                                                                                                                        |

|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|-----|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 6.  | <p>Tim [102x18] MQ Gradients</p> <p>Tim [102x32] MQ Gradients</p> <p>MQ Gradient System with noise reduction features.</p> <p>Integrated measures comprehensively reduce acoustic noise without compromising strongest gradient performance.</p>                                                                                                                                                                                                                                                                                                                                                                                          |
| 7.  | <p>Tim [102x18] MQ Gradients</p> <p>Tim [102x32] MQ Gradients</p> <p>Total imaging matrix with 102 channels (coil elements) that can be connected simultaneously, as well as 18 (32) independent RF channels that can be used simultaneously in one single scan and in one single FOV, each generating an independent partial image. It is for demanding high-end applications and optimized throughput. Tim [102x18] has flexibility in Parallel Imaging. Maximum SNR is ensured through the matrix coil and mMR coil technology.</p>                                                                                                    |
| 8.  | <p>mMR Tim Table</p> <p>Free-floating for PET attenuation optimized table with max. scan range of 200 cm and max. patient weight with vertical movement of 200 kg.</p> <p>Free-floating for PET attenuation optimized table with max. scan range of 200 cm and max. patient weight with vertical movement of 200 kg.</p> <p>Two ergonomically designed and positioned Tableside Control Units, integrated left and right into the front covers.</p> <p>The cantilevered table design gives the system a light and unimposing appearance while providing unobstructed foot space for attending staff and better access to the patient.</p> |
| 9.  | <p>Electronics cabinets:</p> <ul style="list-style-type: none"> <li>- Gradient cabinet</li> <li>- Control cabinet</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| 10. | <p>Separator #3T</p> <p>The SEP (Separation cabinet) has to be used if a central hospital chilled water supply is available or if a chiller of any brand/type is already available.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| 11. | <p>Alarm box</p> <p>The alarm box has the following functions:</p> <ul style="list-style-type: none"> <li>- displays alarm signals</li> <li>- switches the MR system on and off</li> <li>- magnet Stop/ Magnet Quench</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                          |
| 12. | <p>Acquisition Workplace</p> <p>syngo Acquisition Workplace is the workplace in the control room. It includes the host processor with the operating elements monitor and mouse.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| 13. | <p>The image reconstruction computer offers faster processing power for intensive algorithms, amount of data storage for large data sets acquired over long-term measurements, a large amount of main memory for fast processing of measurement data, and a general-purpose graphic processing unit for highly intensive computational calculations.</p>                                                                                                                                                                                                                                                                                  |

|            |                                                                                                                                                                                                                                                                                                                                                                                        |
|------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>14.</b> | <p>Biograph mMR Workflow Suite:</p> <ul style="list-style-type: none"> <li>- Composing including MR-PET composer</li> <li>Inline Composing</li> <li>Inline Diffusion</li> <li>Tim Planning Suite</li> </ul> <p>The Biograph mMR Workflow Suite enables for the first time true simultaneous whole-body MR-PET scanning for head-to-toe imaging.</p>                                    |
| <b>15.</b> | <p>mMR Head/Neck, A Tim Coil</p> <p>The multi-element Matrix coil technology is an essential part supplementing the most innovative Total imaging matrix. Matrix coils have multiple receive coil elements that can be clustered in groups. Each receive coil element is equipped with a low noise preamplifier to maximize signal-to-noise ratio.</p>                                 |
| <b>16.</b> | <p>mMR Head/Neck, A Tim Coil additional (if needed)</p> <p>The multi-element Matrix coil technology is an essential part supplementing the most innovative Total imaging matrix. Matrix coils have multiple receive coil elements that can be clustered in groups. Each receive coil element is equipped with a low noise preamplifier to maximize signal-to-noise ratio.</p>          |
| <b>17.</b> | <p>mMR Spine, A Tim Coil</p> <p>The multi-element Matrix coil technology is an essential part supplementing the most innovative Total imaging matrix. Matrix coils have multiple receive coil elements that can be clustered in groups. Each receive coil element is equipped with a low noise preamplifier to maximize signal-to-noise ratio.</p>                                     |
| <b>18.</b> | <p>mMR Spine, A Tim Coil additional (if needed)</p> <p>The multi-element Matrix coil technology is an essential part supplementing the most innovative Total imaging matrix. Matrix coils have multiple receive coil elements that can be clustered in groups. Each receive coil element is equipped with a low noise preamplifier to maximize signal-to-noise ratio.</p>              |
| <b>19.</b> | <p>mMR Body, A Tim Coil</p> <p>The multi-element PET optimized Matrix coil technology is an essential part supplementing the most innovative Total imaging matrix. Matrix coils have multiple receive coil elements that can be clustered in groups. Each receive coil element is equipped with a low noise preamplifier to maximize signal-to-noise ratio.</p>                        |
| <b>20.</b> | <p>mMR Body, A Tim Coil additional (if needed)</p> <p>The multi-element PET optimized Matrix coil technology is an essential part supplementing the most innovative Total imaging matrix. Matrix coils have multiple receive coil elements that can be clustered in groups. Each receive coil element is equipped with a low noise preamplifier to maximize signal-to-noise ratio.</p> |

|            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>21.</b> | <p>Flex Coil Interface #V,mMR</p> <p>Interface with integrated preamplifiers for the connection of the following flexible coils to the system:</p> <ul style="list-style-type: none"> <li>- 4-channel Flex Coil, large</li> <li>- 4-channel Flex Coil, small</li> <li>- Loop Coil, large</li> <li>- Loop Coil, medium</li> <li>- Loop Coil, small</li> </ul>                                                                                                                                                                                             |
| <b>22.</b> | <p>Flex Coil Interface #V,mMR additional (if needed)</p> <p>Interface with integrated preamplifiers for the connection of the following flexible coils to the system:</p> <ul style="list-style-type: none"> <li>- 4-channel Flex Coil, large</li> <li>- 4-channel Flex Coil, small</li> <li>- Loop Coil, large</li> <li>- Loop Coil, medium</li> <li>- Loop Coil, small</li> </ul>                                                                                                                                                                      |
| <b>23.</b> | <p>Flex Coil, large #V,mMR</p> <p>Lightweight, very flexible, iPAT compatible, 4-element no-tune receiver coil made of soft and smooth material. It is used for examinations of extremities (e.g. medium to large shoulder, hip or knee) and for abdominal examinations.</p> <p>Flex Coil, small #V,mMR</p> <p>Lightweight, very flexible, iPAT compatible, 4-element no-tune receiver coil made of soft and smooth material. It is used for examinations of the upper and lower extremities (e.g. small to medium shoulder, wrist, elbow or ankle).</p> |
| <b>24.</b> | <p>Flex Coil, large #V,mMR(if needed)</p> <p>Lightweight, very flexible, iPAT compatible, 4-element no-tune receiver coil made of soft and smooth material. It is used for examinations of extremities (e.g. medium to large shoulder, hip or knee) and for abdominal examinations.</p>                                                                                                                                                                                                                                                                  |
| <b>25.</b> | <p>Flex Coil, small #V,mMR(if needed)</p> <p>Lightweight, very flexible, iPAT compatible, 4-element no-tune receiver coil made of soft and smooth material. It is used for examinations of the upper and lower extremities (e.g. small to medium shoulder, wrist, elbow or ankle).</p>                                                                                                                                                                                                                                                                   |
| <b>26.</b> | <p>mMR Head 32, A 3T Tim Coil</p> <p>The 32-channel Head Array Coil is an iPAT compatible "no-tune" receiver coil for the Biograph mMR. It can be used for simultaneous MR-PET imaging as well as MR only imaging. A mu-Map for PET attenuation correction of the coil is provided.</p>                                                                                                                                                                                                                                                                  |

|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|-----|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 27. | <p>mMR Breast Coil (if needed)</p> <p>The mMR Breast is used for breast imaging. It consists of two parts that are positioned on and into the frame. Depending on the requirements, the coil may be used in the Head First or Feet First orientation. However, for examinations in the Feet First orientation, the "Whole Body" option is required.</p> <p>The mMR Breast can be used in MR-PET mode, but only in Head First orientation.</p>                                                                                                                            |
| 28. | <p>Tx/Rx CP Extremity Coil (if needed)</p> <p>The TxRx CP Extremity is used for examining knees.</p> <p>Transmit/receive coil for joint examinations in the region of the lower extremities.</p>                                                                                                                                                                                                                                                                                                                                                                         |
| 29. | <p>Small Shoulder Array Coil (if needed)</p> <p>The Shoulder Array is used for examining the shoulder.</p> <p>Receive array coil, anatomically adapted coil tops attached to a base plate, either for the left or for the right shoulder.</p> <p>Base plate is included in scope of delivery for the coil.</p>                                                                                                                                                                                                                                                           |
| 30. | <p>Large Shoulder Array Coil (if needed)</p> <p>The Shoulder Array is used for examining the shoulder.</p> <p>Receive array coil, anatomically adapted coil tops attached to a base plate, either for the left or for the right shoulder.</p> <p>Base plate is included in scope of delivery for the coil.</p>                                                                                                                                                                                                                                                           |
| 31. | <p>Physiological Measurement Unit (PMU)</p> <p>Wireless Physio Control for wireless triggering, synchronizes the measurement with the physiological cycles of cardiac and/or respiratory motion. Wireless technology for all sensors allows fast and easy patient set-up and comfort, and robust cardiac or respiratory signal transmission as it eliminates the need to attach cables to the patient. The Wireless Physio Control contains wireless VCG, respiration and pulse sensors and a charging station as all sensors are powered by rechargeable batteries.</p> |
| 32. | <p>PMU Sensor Kit #Tim</p> <p>Additional PERU, PPU and charging stand.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| 33. | <p>Connection cables</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| 34. | <p>Whole-body superconductive 3T magnet with active shielding (AS) technology with counter coils.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| 35. | <p>Additional Set of Manuals</p> <p>Additional set of manuals for the above selected MR system.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |

#### Accessories

|    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|----|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1. | <p>syngo MR Workplace #mMR</p> <p>syngo MR Workplace based on the newest generation syngo host computer consisting of:</p> <ul style="list-style-type: none"> <li>- Color LCD monitor, 19" flat screen.</li> <li>- syngo MR Workplace computer system with Quad Core Xeon™, 2.66 GHz, 6 GB RAM, dual process</li> <li>host computer incl. one 300 GB system hard disk, one 300 GB hard disk for image data, one CD/DV</li> <li>R drive for image storage on CD/DVD-R and electronic mouse.</li> <li>- syngo MR Maestro class user interface for optimized clinical workflow.</li> </ul> |
| 2. | <p>syngo.via L-Software VB30</p> <p>The syngo.via L-Software offers 2D, 3D, 4D multi-modality routine reading capabilities and a variety of advanced applications tailored to the L-Server HW grade. The combination of syngo.via L-Software and L-Server Hardware is ideal for 2 - 7 users. The availability of all applications and workflows included in syngo.via L-SW is virtually unlimited, i.e. the number of opened cases is only constrained by server HW resources.</p>                                                                                                      |
| 3. | <p>syngo.via Workstation Software VB30</p> <p>The syngo.via Workstation offers 2D, 3D, 4D multi-modality routine reading capabilities and a variety of advanced applications tailored to the Workstation. The combination of syngo.via Software and Workstation Hardware is ideal for 1 - 2 users. The availability of all applications and workflows included in syngo.via Workstation is virtually unlimited, i.e. the number of opened cases is only constrained by server HW resources.</p>                                                                                         |
| 4. | <p>syngo.mMR General #1+</p> <p>syngo.mMR General for one additional user.</p> <p>syngo.mMR General provides additional functionalities within the MR Oncology workflow, for evaluating positron emissions tomography data (PET) in conjunction with MR examinations.</p>                                                                                                                                                                                                                                                                                                               |
| 5. | <p>syngo.MR General Engine</p> <p>syngo.MR General Engine extends syngo.via by adding dedicated workflows and tools for routine and advanced reading of MR examinations.</p>                                                                                                                                                                                                                                                                                                                                                                                                            |
| 6. | <p>Network Switch for RT Interface#T+D</p> <p>Contains network switch for using the software interface in Realtime Interface for connecting to a 3 party computer.</p>                                                                                                                                                                                                                                                                                                                                                                                                                  |
| 7. | <p>Monitor EIZO MX232W col. 2.1MP</p> <p>The EIZO MX 232W is a color widescreen LCD monitor for diagnostic use and clinical review with a resolution of 1920 x 1080 pixels.</p>                                                                                                                                                                                                                                                                                                                                                                                                         |
| 8. | <p>Server HW Config L</p> <p>syngo.via server hardware configuration L.</p> <p>Hewlett Packard rack mount server.</p> <p>syngo.via is a software only medical device. It defines minimum requirements to the hardware it runs on. The hardware itself is not seen as a medical device.</p>                                                                                                                                                                                                                                                                                              |
| 9. | <p>HP Rack 14 Units 19 inch</p> <p>HP Rack Type Rittal for syngo.via server configurations.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |

|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-----|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 10. | UPS 1.5kVA Rackmount 100-127V<br>1.5 kVA 1U Rack mountable Uninterruptible Power Supply for Server usage.<br>UPS 1.5kVA Rackmount 220-240V<br>1.5 kVA 1U Rack mountable Uninterruptible Power Supply for Server usage.                                                                                                                                                                                                                                                                                                                                                                  |
| 11. | Monitor for Administration<br>HP LCD monitor 20" for syngo administration.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| 12. | Workplace/Workstation Hardware<br>syngo.via Server-based Workstation HW, tower floorstand configuration.<br>This server solution is an entry level solution. The maximum number of concurrent users for advanced applications* on the workstation HW is limited to one user (see syngo.via Data Sheet).<br>The syngo.via Workstation SW is ideal for 1-2 users with 16k max. # of slices for concurrent rendering. syngo.via CT/MI/RT and BreastCare Workplace SW is limited to 1 user.<br>For a higher redundancy, please order the Server HW config L or larger server configuration. |
| 13. | Ultramini Workstation Client<br>syngo.via Client Hardware based on small form factor workstation from HP Inc                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| 14. | syngo.via BreastCare WP SW VB30<br>syngo.via BreastCare Workplace is a dedicated reading environment for mammography (DICOM MG). Together with high-resolution monitors approved for mammography this workplace provides excellent image quality and efficient workflow for breast cancer diagnostic and follow-up. syngo.via BreastCare Workplace is multi-vendor compatible according the IHE profile Mammography Image Display.<br>syngo.via BreastCare Workplace is available as a one user workplace configuration.                                                                |
| 15. | The Biograph mMR Workflow Suite enables for the first time true simultaneous whole-body MR-PET scanning for head-to-toe imaging.                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| 16. | Angio Dot Engine #mMR<br>The timing of contrast injection and scan is widely considered the most challenging part of an angiographic exam. Angio Dot guides the user through angiographic single or multi station examinations by providing semi-automatic detection of arterial and venous timing windows using a test bolus technique.                                                                                                                                                                                                                                                |
| 17. | TimCT Angio Dot Engine #T+D<br>syngo TimCT Angiography employs the revolutionary TimCT Continuous Table move technology for large Field of View angiographies with the smoothest workflow and the most homogeneous image quality. syngo TimCT Angiography is built on the Tim technology as well as on a highly advanced patient table with high positioning accuracy and an RF shielded table drive.                                                                                                                                                                                   |
| 18. | TimCT Angio Dot Engine #T+D<br>syngo TimCT Angiography employs the revolutionary TimCT Continuous Table move technology for large Field of View angiographies with the smoothest workflow and the most homogeneous image quality. syngo TimCT Angiography is built on the Tim technology as well as on a highly advanced patient table with high positioning accuracy and an RF shielded table drive.                                                                                                                                                                                   |
| 19. | TimCT Onco Dot Engine #mMR<br>TimCT Onco Dot engine adds the excellent Dot workflow functionality to TimCT Oncology features. (The TimCT Oncology features are standard with Biograph E11P).                                                                                                                                                                                                                                                                                                                                                                                            |

|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|-----|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 20. | <p>TimCT Onco Dot Engine #mMR</p> <p>TimCT Onco Dot engine adds the excellent Dot workflow functionality to TimCT Oncology feature (The TimCT Oncology features are standard with Biograph E11P).</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| 21. | <p>Abdomen Dot Engine #T+D</p> <p>The Abdomen Dot Engine offers a comprehensive set of guidance and automation, so that robust image quality can be achieved fast and independently from the user.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| 22. | <p>Prostate Package #T+D</p> <p>The prostate spectroscopy package is a comprehensive software package which bundles:</p> <ul style="list-style-type: none"> <li>- Single Voxel Spectroscopy</li> <li>- 2D Chemical Shift Imaging</li> <li>- 3D Chemical Shift Imaging</li> <li>- Spectroscopy Evaluation syngo</li> <li>- syngo Tissue 4D Evaluation</li> </ul>                                                                                                                                                                                                                                                                                                                            |
| 23. | <p>Tissue 4D syngo #Tim</p> <p>Tissue 4D is an application for visualizing and post-processing dynamic contrast-enhanced 3D datasets.</p> <p>This card provides two evaluation options:</p> <ul style="list-style-type: none"> <li>- Standard curve evaluation</li> <li>- Curve evaluation according to a pharmacokinetic model</li> </ul> <p>Tissue 4D syngo MRWP #Tim</p> <p>Tissue 4D is an application for visualizing and post-processing dynamic contrast-enhanced 3D datasets.</p> <p>This card provides two evaluation options:</p> <ul style="list-style-type: none"> <li>- Standard curve evaluation</li> <li>- Curve evaluation according to a pharmacokinetic model</li> </ul> |
| 24. | <p>syngo.MR Onco Engine #1</p> <p>syngo.MR Onco Engine pools the functions from syngo.MR Onco and from syngo.MR 3D Lesion Segmentation for efficient oncology evaluation and analysis.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| 25. | <p>syngo.MR Onco Engine #1+</p> <p>License for one additional user for syngo.MR Onco Engine.</p> <p>syngo.MR Onco Engine pools the functions from syngo.MR Onco and from syngo.MR 3D Lesion Segmentation for efficient oncology evaluation and analysis.</p>                                                                                                                                                                                                                                                                                                                                                                                                                               |
| 26. | <p>syngo.MR Onco Engine Pro #1</p> <p>syngo.MR Onco Engine Pro, in conjunction with syngo.MR Onco Engine, provides comprehensive functionalities for the analysis of oncological MR studies. syngo.MR Onco Engine Pro includes the functionalities of syngo.MR OncoCare as an enhancement package for syngo.MR Onco Engine.</p>                                                                                                                                                                                                                                                                                                                                                            |
| 27. | <p>syngo.MR Onco Engine Pro #1+</p> <p>License for one additional user for syngo.MR Onco Engine Pro.</p> <p>syngo.MR Onco Engine Pro, in conjunction with syngo.MR Onco Engine, provides comprehensive functionalities for the analysis of oncological MR studies. syngo.MR Onco Engine Pro includes the functionalities of syngo.MR OncoCare as an enhancement package for syngo.MR Onco Engine.</p>                                                                                                                                                                                                                                                                                      |

|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|-----|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 28. | <p><b>TWIST syngo #Tim</b></p> <p>This package contains a Siemens unique sequence and protocols for time-resolved (4D) MR angiographic and dynamic imaging in general with high spatial and temporal resolution. syngo TWIST supports comprehensive dynamic MR angio exams in all body regions. It offers temporal information of vessel filling in addition to conventional static MR angiography, which can be beneficial in detecting or evaluating malformations such as shunts.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| 29. | <p><b>syngo.MR Vascular Analysis #1</b></p> <p>syngo.MR Vascular Analysis for fast and intuitive assessment/ quantification of general vascular pathologies, such as stenosis.</p> <p><b>Vessel View syngo #Tim</b></p> <p>syngo Vessel View provides a dedicated evaluation software for automated image analysis (e.g. stenosis quantification) of MR and CT angiography data for small and large vessels.</p> <p><b>Vessel View AVS #Tim</b></p> <p>As a separate additional option of Vessel View, this package allows for automated segmentation and separation of arteries and veins, as well as suppression of surrounding tissue. Supports modes allowing the display of only arteries or only veins, or arteries and veins together in different colors.</p> <p><b>Vessel View syngo MRWP #Tim</b></p> <p>This application provides a dedicated evaluation software for automated image analysis (e.g. stenosis quantification) of MR and CT angiography data for small and large vessels.</p> <p><b>Vessel View AVS MRWP #Tim</b></p> <p>As a separate additional option of Vessel View, this package allows for automated segmentation and separation of arteries and veins, as well as suppression of surrounding tissue. Supports modes allowing the display of only arteries or only veins, or arteries and veins together in different colors.</p> |
| 30. | <p><b>TWIST syngo #Tim</b></p> <p>This package contains a Siemens unique sequence and protocols for time-resolved (4D) MR angiographic and dynamic imaging in general with high spatial and temporal resolution. syngo TWIST supports comprehensive dynamic MR angio exams in all body regions. It offers temporal information of vessel filling in addition to conventional static MR angiography, which can be beneficial in detecting or evaluating malformations such as shunts.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| 31. | <p><b>BLADE</b></p> <p>Motion and flow insensitive multi-shot Turbo Spin Echo (TSE) sequence for all body regions with optional inter-shot motion correction.</p> <p>BLADE supports T2-weighted, dark fluid and STIR contrast imaging as well as inversion recovery T1-weighted imaging.</p> <p>BLADE uses an acquisition scheme which is less sensitive to motion and flow than the scheme used by conventional TSE. In the head additionally a 2D (in-plane) motion correction can be used to further reduce artifacts resulting from patient motion during the scan.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| 32. | <p><b>Cardiac Dot Engine #T+D</b></p> <p>Cardiac examinations used to be the most complex exams in MR. Now Cardiac Dot Engine supports the user in many ways. Using anatomical landmarks, standard views of the heart, such as dedicated long axis and short-axis views, are easily generated and can easily be reproduced using different scanning techniques.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| 33. | <p><b>syngo.MR Cardio Engine #1</b></p> <p>The syngo.MR Cardio Engine bundles the following features for Cardiac evaluation:</p> <ul style="list-style-type: none"> <li>- syngo.MR Cardiac 4D Ventricular Function</li> <li>- syngo.MR Cardiac Flow</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |

|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|-----|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 34. | <p>syngo.MR Cardio Engine #1+</p> <p>License for one additional user for syngo.MR Cardio Engine.</p> <p>The syngo.MR Cardio Engine bundles the following features for Cardiac evaluation:</p> <ul style="list-style-type: none"> <li>- syngo.MR Cardiac 4D Ventricular Function</li> <li>- syngo.MR Cardiac Flow</li> </ul>                                                                                                                                                                                                                                                                                                                    |
| 35. | <p>Advanced Cardiac Package #T+D</p> <p>This package contains special sequences and protocols for advanced cardiac imaging including 3D 4D syngo BEAT functionalities. It supports advanced techniques for ventricular function imaging, dynamic imaging, tissue characterization, coronary imaging, and more.</p> <p>MyoMaps #3T</p> <p>This package contains special sequences and protocols for inline T1 and T2 calculation at the heart. The generation of T1 and T2 parametric maps is enhanced by the use of motion correction. T1 and T2 parametric maps could be used to support assessment of cardiovascular disease.</p>            |
| 36. | <p>HD Chest #mMR</p> <p>Expands the Biograph mMR PET Respiratory option to include HD-Chest technology. HD-Chest uses adaptive respiratory gating to reduce motion artifacts while providing optimal count statistics.</p> <p>HD Chest MRWP #mMR</p> <p>Expands the Biograph mMR PET Respiratory option to include HD-Chest technology. HD-Chest uses adaptive respiratory gating to reduce motion artifacts while providing optimal count statistics.</p>                                                                                                                                                                                     |
| 37. | <p>Multi-gate #mMR</p> <p>Expands the Biograph mMR PET Respiratory option to provide multiple PET respiratory gates, with matching, simultaneously acquired MR gates for attenuation correction.</p> <p>Multi-gate MRWP #mMR</p> <p>Expands the Biograph mMR PET Respiratory option to provide multiple PET respiratory gates, with matching, simultaneously acquired MR gates for attenuation correction.</p>                                                                                                                                                                                                                                 |
| 38. | <p>Flow Quantification #Tim</p> <p>Special sequences for quantitative assessment of flow.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| 39. | <p>Argus Flow</p> <p>Evaluation software for automated image analysis of blood and cerebrospinal fluid (CSF) flow.</p> <p>Argus Flow MRWP</p> <p>Evaluation software for automated image analysis of blood and cerebrospinal fluid (CSF) flow.</p>                                                                                                                                                                                                                                                                                                                                                                                             |
| 40. | <p>NATIVE syngo</p> <p>Integrated software package with sequences and protocols for non-contrast enhanced 3D MRA with high spatial resolution. syngo NATIVE particularly enables imaging of abdominal and peripheral vessels and is an alternative to MR angiography techniques with contrast medium, especially for patients with severe renal insufficiency.</p>                                                                                                                                                                                                                                                                             |
| 41. | <p>Argus DynamicSignal</p> <p>Evaluation software for automated image analysis of dynamic studies in the body. The combination of automated contouring, motion correction and easy-to-use editing tools, provides users with a rapid way to quantify regional changes in the tissue based on signal-time curves.</p> <p>Argus DynamicSignal MRWP</p> <p>Evaluation software for automated image analysis of dynamic studies in the body. The combination of automated contouring, motion correction and easy-to-use editing tools, provides users with a rapid way to quantify regional changes in the tissue based on signal-time curves.</p> |

|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|-----|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 42. | Argus 4D Ventr. Function syngo #Tim<br>syngo Argus 4D Ventricular Function software processes MR cine images of the heart and generates quantitative results for physicians in the diagnostic process.<br>Argus 4D Ventr.Funct.syngo MRWP#Tim<br>syngo Argus 4D Ventricular Function software processes MR cine images of the heart and generates quantitative results for physicians in the diagnostic process.                                                                                                                            |
| 43. | Dynamic Analysis MR evaluation software allows the calculation of functions such as addition/subtraction, division/multiplication, ADC maps, T1 and T2, z-Score (t-Test), Time-to-Peak maps (TTP) and standard deviation                                                                                                                                                                                                                                                                                                                    |
| 44. | Argus viewer can be used to display MR cine studies. Argus Viewer allows users to load a large list of dynamic data sets and view them comfortably. This is a feature that greatly reduces the reading and review time for MR cardiac studies.                                                                                                                                                                                                                                                                                              |
| 45. | Interactive Realtime #T+D<br>Interactive Realtime Imaging for interactive acquisition, reconstruction and display of image data. Uses ultra-fast TrueFISP and other gradient-echo sequences of high image contrast. The user can interactively navigate in all planes on-the-fly during data acquisition.<br>Realtime Interface #T+D<br>Adds advanced features to Interactive Realtime, including the capability for multi-slice acquisition, definition of acquisition order, pausing, mosaic display, and skipping of the physio trigger. |
| 46. | CAIPIRINHA<br>Enables VIBE sequence with improved iPAT <sup>2</sup> algorithm to improve abdominal dynamic scans as well as SNR.<br>Reduced patient stress can be achieved through reduced acquisition (and breathhold) times.                                                                                                                                                                                                                                                                                                              |
| 47. | 1D/2D PACE (Prospective Acquisition CorrEction) – the motion freeze for examinations with breath-hold and free breathing combined with PET.                                                                                                                                                                                                                                                                                                                                                                                                 |
| 48. | Brain Compass #mMR<br>This feature is the first Motion Correction application acquired simultaneously. Brain Compass combines functional MRI with Motion corrected static PET.<br>Brain Compass MRWP #mMR<br>This feature is the first Motion Correction application acquired simultaneously. Brain Compass combines functional MRI with Motion corrected static PET.                                                                                                                                                                       |
| 49. | The Spine Dot Engine provides optimized cervical, thoracic and lumbar spine imaging. Amongst various features to support a streamlined spine workflow is AutoLabeling of the vertebrae. The Spine Dot Engine provides optimized cervical, thoracic and lumbar spine imaging for patients of all conditions.<br>Spine Dot Engine provides the functionality to simplify your spine workflow by providing tools to reduce examination times, achieve optimal image quality, and assist you during reading.                                    |
| 50. | Large Joint Dot Engine E11 #T+D<br>Large Joint Dot Engine optimizes image quality of knee, hip and shoulder scans by proposing the most appropriate protocols according to the examination strategy chosen for the specific patient. It ensures reproducible image quality and streamlines large joint examinations to the greatest extent.                                                                                                                                                                                                 |
| 51. | Composing syngo MRWP #Tim<br>This application provides dedicated evaluation software for creation of full-format images from overlapping MR volume data sets and MIPs acquired at multiple stages.                                                                                                                                                                                                                                                                                                                                          |

|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|-----|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 52. | <p>The Spine Dot Engine provides optimized cervical, thoracic and lumbar spine imaging. Amongst various features to support a streamlined spine workflow is AutoLabeling of the vertebrae. The Spine Dot Engine provides optimized cervical, thoracic and lumbar spine imaging for patients in all conditions.</p> <p>Spine Dot Engine provides the functionality to simplify your spine workflow by providing tools to reduce examination times, achieve optimal image quality, and assist you during reading.</p> <ul style="list-style-type: none"> <li>- User guidance step-by-step</li> </ul> <p>AutoPosition<br/>AutoAlign Spine with intervertebral disc detection<br/>AutoCoverage<br/>AutoSatPosition<br/>Initial and interactive snapping<br/>AutoLabeling of vertebrae</p> <p>The Spine Dot Engine includes:</p> <ul style="list-style-type: none"> <li>- syngo WARP - Susceptibility Artifact Reduction syngo WARP integrates different techniques tailored to reduce susceptibility artifacts caused by orthopedic MR-conditional metal implants. 2D TSE sequence combining optimized high-bandwidth protocols and View Angle Tilting (VAT) techniques. This helps in evaluation of soft tissue in proximity of the implant. Available protocols include T1-, T2-weighted, and STIR contrast.</li> </ul> |
| 53. | <p>The Spine Dot Engine provides optimized cervical, thoracic and lumbar spine imaging. Amongst various features to support a streamlined spine workflow is AutoLabeling of the vertebrae. The Spine Dot Engine provides optimized cervical, thoracic and lumbar spine imaging for patients in all conditions.</p> <p>Spine Dot Engine provides the functionality to simplify your spine workflow by providing tools to reduce examination times, achieve optimal image quality, and assist you during reading.</p> <ul style="list-style-type: none"> <li>- User guidance step-by-step</li> </ul> <p>AutoPosition<br/>AutoAlign Spine with intervertebral disc detection<br/>AutoCoverage<br/>AutoSatPosition<br/>Initial and interactive snapping<br/>AutoLabeling of vertebrae</p> <p>The Spine Dot Engine includes:</p> <ul style="list-style-type: none"> <li>- syngo WARP - Susceptibility Artifact Reduction syngo WARP integrates different techniques tailored to reduce susceptibility artifacts caused by orthopedic MR-conditional metal implants. 2D TSE sequence combining optimized high-bandwidth protocols and View Angle Tilting (VAT) techniques. This helps in evaluation of soft tissue in proximity of the implant. Available protocols include T1-, T2-weighted, and STIR contrast.</li> </ul> |
| 54. | <p>3D distortion correction filter (inline and offline) for high spatial accuracy DICOM MPPS (Modality Performed Procedure Steps) for efficient organization of workflow within HIS/RIS systems. MPPS allows to communicate examination data from the MR system to an information system (e.g. RIS system) and to provide data for billing, documentation and planning purposes to the information system.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| 55. | <p>Soft Tissue Motion Correction #Tim<br/>3D post-processing step, elastic soft tissue/breast motion correction.<br/>Soft Tissue Motion Corr. MRWP #Tim<br/>3D post-processing step, elastic soft tissue/breast motion correction.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |

|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|-----|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 56. | <p>syngo THERM</p> <p>syngo THERM supports the operator to identify areas that show a high amount of temperature change in non-moving organs.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| 57. | <p>Advanced Diffusion #mMR</p> <p>QuietX DWI enables quieter diffusion-weighted imaging of the brain with up to 70% reduction in sound pressure relative to conventional diffusion-weighted imaging.</p> <p>RESOLVE (Readout Segmentation Of Long Variable Echo-trains) is a multi-shot, readout segmented EPI sequence for high-resolution, low-distortion diffusion-weighted imaging (DWI). This technique is largely insensitive to susceptibility effects, providing anatomically accurate diffusion imaging for the brain, spine, breast and prostate. In combination with the DTI Tractography package, RESOLVE enables excellent</p>                                                                                                                                                                                                                                                                                                                                                                        |
| 58. | <p>Diffusion Tensor Imaging #mMR</p> <p>Diffusion Tensor Imaging provides a Single Shot EPI sequence for measuring diffusion-weighted data sets with up to 256 directions of diffusion weighting. Based on these data sets, the diffusion tensor itself and parametric maps derived from it (e.g. fractional anisotropy) are calculated automatically and in real-time.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| 59. | <p>DTI Evaluation #Tim</p> <p>DTI Evaluation provides advanced post-processing and visualization of Diffusion Tensor Imaging (DTI) data. DTI Evaluation includes the possibility to calculate tensor data from a DTI dataset retrospectively and enables calculation of different diffusion parameter maps. Furthermore, it facilitates joint ROI-based evaluation of parameter images and anatomical images, as well as color-coded display and fused 3D visualization in the anatomical context.</p> <p>DTI Evaluation MRWP #Tim</p> <p>DTI Evaluation provides advanced post-processing and visualization of Diffusion Tensor Imaging (DTI) data. DTI Evaluation includes the possibility to calculate tensor data from a DTI dataset retrospectively and enables calculation of different diffusion parameter maps. Furthermore, it facilitates joint ROI-based evaluation of parameter images and anatomical images, as well as color-coded display and fused 3D visualization in the anatomical context.</p> |
| 60. | <p>DTI Tractography syngo #Tim</p> <p>syngo DTI Tractography allows the visualization of multiple white matter tracts of the human brain based on diffusion tensor imaging data. DTI Tractography is optimized to support the presurgical planning and to allow for neuro physiological research with respect to connectivity and white matter pathology.</p> <p>DTI Tractography syngo MRWP #Tim</p> <p>syngo DTI Tractography allows the visualization of multiple white matter tracts of the human brain based on diffusion tensor imaging data. DTI Tractography is optimized to support the presurgical planning and to allow for neuro physiological research with respect to connectivity and white matter pathology.</p>                                                                                                                                                                                                                                                                                   |
| 61. | <p>DTI Package #mMR</p> <p>The DTI Package is a bundle of:</p> <ul style="list-style-type: none"> <li>- Diffusion Tensor Imaging</li> <li>- DTI Evaluation and</li> <li>- DTI Tractography syngo</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |

|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|-----|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 62. | <p><b>DTI Package #mMR</b></p> <p>The DTI Package is a bundle of:</p> <ul style="list-style-type: none"> <li>- Diffusion Tensor Imaging</li> <li>- DTI Evaluation and</li> <li>- DTI Tractography syngo</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| 63. | <p><b>Neuro Perfusion Package #T+D</b></p> <p>The Neuro Perfusions Package helps to streamline the clinical workflow by inline post-processing of dynamic susceptibility contrast (DSC) based perfusion imaging. This makes it possible to see perfusion maps immediately.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| 64. | <p><b>Neuro Perfusion Package #T+D</b></p> <p>The Neuro Perfusions Package helps to streamline the clinical workflow by inline post-processing of dynamic susceptibility contrast (DSC) based perfusion imaging. This makes it possible to see perfusion maps immediately.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| 65. | <p><b>Neuro Perfusion Evaluation #T+D</b></p> <p>Neuro Perfusion Evaluation syngo provides a task card for detailed post-processing of brain perfusion data sets. Color display of the relative Mean Transit Time (relMTT), relative Cerebral Blood Volume (relCBV), corrected rel CBV, and relative Cerebral Blood Flow (relCBF) is supported. Flexible selection of the Arterial Input Function (AIF) for more reliable analysis taking into account the dynamics over time of the contrast agent enhancement. Furthermore, a calculation of maps using automatically selected local Arterial Input Functions (AIF) is provided to reduce the amount of user interactions.</p> <p><b>Neuro Perfusion Evaluation MRWP#T+D</b></p> <p>Neuro Perfusion Evaluation syngo provides a task card for detailed post-processing of brain perfusion data sets. Color display of the relative Mean Transit Time (relMTT), relative Cerebral Blood Volume (relCBV), corrected rel CBV, and relative Cerebral Blood Flow (relCBF) is supported. Flexible selection of the Arterial Input Function (AIF) for more reliable analysis taking into account the dynamics over time of the contrast agent enhancement. Furthermore a calculation of maps using automatically selected local Arterial Input Functions (AIF) is provided to reduce the amount of user interactions.</p> |
| 66. | <p><b>Arterial Spin Labeling 2D</b></p> <p>ASL is a non-contrast enhanced brain perfusion technique. EPI sequence enhanced for ASL (Arterial Spin Labeling) with preparation module (inversion pulse, saturation pulses) and selectable prospective motion correction. Perfusion-weighted color maps and relative cerebral blood flow (relCBF) color maps are calculated with Inline technology.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| 67. | <p><b>Arterial Spin Labeling 3D #T+D</b></p> <p>ASL is a non-contrast enhanced brain perfusion technique. A 3D volume is acquired with high SNR using a turbo gradient spin echo technique and an ASL preparation module to achieve clinically feasible scan times.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| 68. | <p><b>SWI #Tim</b></p> <p>Susceptibility Weighted Imaging is a high-resolution 3D imaging technique for the brain with ultra-high sensitivity for microscopic magnetic field inhomogeneities caused by deoxygenated blood, products of blood decomposition and microscopic iron deposits. Among other things, the method allows for the highly sensitive proof of cerebral hemorrhages and the high-resolution display of venous cerebral blood vessels.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| 69. | <p><b>Inline BOLD Imaging #Tim</b></p> <p>The BOLD imaging package is based on blood oxygen level dependent (BOLD) contrast-sensitive single-shot EPI sequences.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |

|            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>70.</b> | <p><b>BOLD 3D Evaluation syngo #Tim</b><br/>BOLD 3D Evaluation syngo is the comprehensive processing and visualization package for BOLD fMRI. It provides a full set of features for clinical fMRI, as well as advanced features for more research-oriented applications.</p> <p><b>BOLD 3D Evaluation syngo MRWP #Tim</b><br/>BOLD 3D Evaluation syngo is the comprehensive processing and visualization package for BOLD fMRI. It provides a full set of features for clinical fMRI, as well as advanced features for more research-oriented applications.</p> |
| <b>71.</b> | <p><b>3D PACE syngo #Tim</b><br/>3D PACE (Prospective Acquisition CorrEction) enhances Inline BOLD imaging with motion correction during the acquisition of a BOLD exam. In contrast to a retrospective motion correction that corrects previously acquired data, the unique 3D PACE tracks the head of the patient, correcting for motion in real time during the acquisition.</p>                                                                                                                                                                              |
| <b>72.</b> | <p><b>Neuro fMRI Package #T+D</b><br/>The Neuro fMRI Package is a bundle of:</p> <ul style="list-style-type: none"> <li>- Inline BOLD Imaging</li> <li>- 3D PACE syngo</li> <li>- BOLD 3D Evaluation syngo</li> </ul>                                                                                                                                                                                                                                                                                                                                            |
| <b>73.</b> | <p><b>Neuro fMRI/DTI Combi Package #mMR</b><br/>The Neuro fMRI/DTI Combi Package is a bundle of:</p> <ul style="list-style-type: none"> <li>- Inline BOLD Imaging</li> <li>- 3D PACE syngo</li> <li>- BOLD 3D Evaluation syngo</li> <li>- Diffusion Tensor Imaging</li> <li>- DTI Evaluation</li> <li>- DTI Tractography syngo</li> </ul>                                                                                                                                                                                                                        |
| <b>74.</b> | <p><b>Spectroscopy Package #T+D</b><br/>The Spectroscopy Package is a comprehensive software package which bundles Single Voxel Spectroscopy, 2D Chemical Shift Imaging, 3D Chemical Shift Imaging and syngo Spectroscopy Evaluation.</p>                                                                                                                                                                                                                                                                                                                        |
| <b>75.</b> | <p><b>Single Voxel Spectroscopy #Tim</b><br/>Integrated software package including sequences and protocols for proton spectroscopy to examine metabolic changes in the brain (e.g. in tumors and degenerative diseases).</p>                                                                                                                                                                                                                                                                                                                                     |
| <b>76.</b> | <p><b>Spectroscopy Evaluation syngo</b><br/>Comprehensive software package for evaluation of spectroscopy data on the syngo Acquisition Workplace.</p> <p><b>Spectroscopy Eval. syngo MRWP</b><br/>Comprehensive software package for evaluation of spectroscopy data on the syngo MR Workplace.</p>                                                                                                                                                                                                                                                             |
| <b>77.</b> | <p><b>2D Chemical Shift Imaging #Tim</b><br/>Sequences and protocols for proton 2D Chemical Shift Imaging (2D CSI).</p>                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>78.</b> | <p><b>3D Chemical Shift Imaging #Tim</b><br/>Sequences and protocols for proton 3D Chemical Shift Imaging (3D CSI).</p>                                                                                                                                                                                                                                                                                                                                                                                                                                          |

|     |                                                                                                                                                                                                                                                                                                                                                                                           |
|-----|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 79. | GRACE syngo #3T<br>Integrated software package including sequences and protocols for proton spectroscopy, optimized for breast studies.                                                                                                                                                                                                                                                   |
| 80. | syngo.MR Neuro Perfusion #1<br>syngo.MR Neuro Perfusion extends the MR Neurology workflow with advanced processing tools for the analysis of brain perfusion datasets                                                                                                                                                                                                                     |
| 81. | syngo.MR Neuro Perfusion #1+<br>License for one additional user for syngo.MR Neuro Perfusion.<br>Contains MR Neuro Perfusion evaluation and Automatic local AIF calculation                                                                                                                                                                                                               |
| 82. | syngo.MR Neuro Perfusion Engine #1<br>syngo.MR Neuro Perfusion Engine extends the MR Neurology workflow with a complete package for advanced processing and evaluation of brain perfusion datasets                                                                                                                                                                                        |
| 83. | syngo.MR Neuro Perfusion Engine #1+<br>License for one additional user for syngo.MR Neuro Perfusion Engine.                                                                                                                                                                                                                                                                               |
| 84. | syngo.MR Neuro 3D Engine #1<br>The package syngo.MR Neuro 3D Engine consists of:<br>syngo.MR Tractography, enabling the representation of diffusion paths of the human brain based on diffusion tensor imaging. syngo.MR Tractography supports operation planning and is suitable for neurophysiological research in relation to cortical networking and pathologies of the white matter  |
| 85. | syngo.MR Neuro 3D Engine #1+<br>The package syngo.MR Neuro 3D Engine consists of:<br>syngo.MR Tractography, enabling the representation of diffusion paths of the human brain based on diffusion tensor imaging. syngo.MR Tractography supports operation planning and is suitable for neurophysiological research in relation to cortical networking and pathologies of the white matter |
| 86. | Multinuclear Support #mMR<br>Hardware and software package to prepare the MR system for spectroscopy and nuclear imaging with the nuclei $^3\text{He}$ , $^7\text{Li}$ , $^{13}\text{C}$ , $^{17}\text{O}$ , $^{19}\text{F}$ , $^{23}\text{Na}$ , $^{31}\text{P}$ , and $^{129}\text{Xe}$ .                                                                                               |
| 87. | Multinuclear Spectroscopy #mMR<br>Software package for CSI spectroscopy with $^7\text{Li}$ , $^{13}\text{C}$ , $^{17}\text{O}$ , $^{19}\text{F}$ , and $^{31}\text{P}$ nuclei.                                                                                                                                                                                                            |
| 88. | 3D VRT syngo<br>This application provides a dedicated evaluation software for colored, volume rendered (VRT) visualization of MR and CT data as a complement to traditional MIP display.<br>3D VRT syngo MRWP<br>This application provides a dedicated evaluation software for colored, volume rendered (VRT) visualization of MR and CT data as a complement to traditional MIP display. |

|                                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|----------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| oscopy, opti                                       | 89. Fly Through syngo<br>Integrated syngo image analysis is a vital component of the MR examination.<br>This application provides dedicated evaluation software for simulation of virtual endoscopy or bronchoscopy and for Fly-Through in vascular structures, where real endoscopic procedures are impossible.<br>Fly Through syngo MRWP<br>Integrated syngo image analysis is a vital component of the MR examination.<br>This application provides dedicated evaluation software for simulation of virtual endoscopy or bronchoscopy and for Fly-Through in vascular structures, where real endoscopic procedures are impossible. |
| rocessing too                                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| complete packa                                     | 90. syngo Image Fusion<br>This application provides a dedicated evaluation software for spatial alignment (matching) and visualization of image data either from different modalities (CT,MR,NM,PET) or from the same modality but from multiple examinations of the same patient. It supports optimal diagnostic outcome (fusion of morphological and functional information) and therapy planning.                                                                                                                                                                                                                                  |
| nan brain bas<br>d is suitable fo<br>the white mat | 91. 3D DESS for high-resolution 3D studies of joints with excellent T1 and T2 contrast. DESS is a double echo gradient sequence in which a FISP (transversally rephased gradient echo) and a PSIF (RF refocused gradient echo) echo are simultaneously acquired and added, resulting in an increased signal-to-noise ratio.                                                                                                                                                                                                                                                                                                           |
|                                                    | 92. 3D CISS for excellent very high-resolution studies, especially useful in inner-ear examinations.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| nan brain bas<br>d is suitable fo<br>the white mat | 93. MapIt syngo #Tim<br>Based on the T1, T2 or T2* properties of the cartilage syngo ParametricMap allows the early detection of osteoarthritic break down of cartilage structures even before morphological changes occur. The method supports therapeutic decisions in individual patients and can be used to control treatments non-invasively, replacing surgeries or biopsies.                                                                                                                                                                                                                                                   |
| nuclear imagi                                      | 94. GMR (Gradient Motion Rephasing). Sequences with additional bipolar gradient pulses, permitting effective reduction of flow artifacts                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|                                                    | 95. Quiet Suite #T+D<br>Quiet Suite enables complete, quiet examinations for neurology and orthopedics with at least 70% reduction in sound pressure levels.                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| ndered (VRT)                                       | 96. 3D distortion correction filter (inline and offline) for high spatial accuracy DICOM MPPS (Modality Performed Procedure Steps) for efficient organization of workflow within HIS/RIS systems. MPPS allows to communicate examination data from the MR system to an information system (e.g. RIS system) and to provide data for billing, documentation and planning purposes to the information system.                                                                                                                                                                                                                           |
| ndered (VRT)                                       | 97. syngo Expert-i #T+D<br>This software application enables remote access to the system (connected via local area network) for planning and processing.<br>syngo Expert-i MRWP #T+D<br>This software application enables remote access to the system (connected via local area network) for planning and processing.                                                                                                                                                                                                                                                                                                                 |

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 98.  | <p>syngo Security Package</p> <p>Software option for general regulatory security rules, providing enhanced security features including user management and audit trail functionality. This package supports customers in their achieving compliance with the HIPAA "Privacy" rule.</p> <p>syngo Security Package MRWP</p> <p>Software option for general regulatory security rules, providing enhanced security features including user management and audit trail functionality. This package supports customers in their achieving compliance with the HIPAA "Privacy" rule.</p> |
| 99.  | <p>Dual Monitor Package #mMR</p> <p>The Dual Monitor is an extended feature which facilitates the display of the e.g. PET Acquisition Monitor for a better overview during the mMR examination.</p>                                                                                                                                                                                                                                                                                                                                                                                |
| 100. | <p>syngo.MR General Engine</p> <p>syngo.MR General Engine extends syngo.via by adding dedicated workflows and tools for routine and advanced reading of MR examinations.</p>                                                                                                                                                                                                                                                                                                                                                                                                       |
| 101. | <p>The syngo.via CTWP Starter Engine includes the most frequently used applications in routine</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| 102. | <p>LiverLab #T+D</p> <p>LiverLab is a system guided workflow to examine the hepatic fat and iron status, as part of the Abdomen Dot Engine.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| 103. | <p>syngo.CT Colonography #1</p> <p>syngo.CT Colonography is a non-invasive, diagnostic tool for CT data that can be used to locate and evaluate lesions in the colon.</p>                                                                                                                                                                                                                                                                                                                                                                                                          |
| 104. | <p>syngo.CT CaScoring #1</p> <p>syngo.CT CaScoring for quick risk assessment and coronary age calculation</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| 105. | <p>syngo.CT Bone Reading #1</p> <p>syngo.CT Bone Reading provides a specific visualization for diagnostic bone assessment while reading radiological CT datasets. All ribs and the spine are displayed spread out in one plane and interactive functionality provides scrolling and rotating. Additionally, labels and numbers attached to the ribs and the spine. The supplemental second reader tool provides CAD (Computer Aided Detection) functionality, which is designed to assist radiologists in the detection of suspicious lesions of the spine.</p>                    |
| 106. | <p>syngo.CT Pulmo 3D #1</p> <p>syngo.via application for CT-based clinical assessment of lung diseases like COPD (Chronic Obstructive Pulmonary Disease) and associated lung emphysema. Provides automated evaluation and documentation by 3D quantification of left and right lung, lung lobes, and automated segmentation and measurements of the trachea and associated bronchi.</p>                                                                                                                                                                                            |
| 107. | <p>syngo.MR Cardiac Perfusion #1</p> <p>syngo.MR Cardiac Perfusion processes cardiac MR perfusion images to evaluate ischemic heart disease e.g. perfusion at rest and perfusion at stress. The application generates pixel based semi-quantitative parametric maps to support the diagnostic process. The MR cardiac interactive report template is included.</p>                                                                                                                                                                                                                 |

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 108. | 3D distortion correction filter (inline and offline) for high spatial accuracy DICOM MPPS (Modality Performed Procedure Steps) for efficient organization of workflow within HIS/RIS systems. MPPS allows to communicate examination data from the MR system to an information system (e.g. RIS system) and to provide data for billing, documentation and planning purposes to the information system.                                                                                                                                                                                                    |
| 109. | syngo.MR General Engine<br>syngo.MR General Engine extends syngo.via by adding dedicated workflows and tools for routine and advanced reading of MR examinations.                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| 110. | syngo.mMR General #1<br>syngo.mMR General provides additional functionalities within the MR Oncology workflow, for evaluating positron emissions tomography data (PET) in conjunction with MR examinations. Results from MR-PET examinations are loaded into suitable layouts and optimally structured. syngo.mMR General allows for quantitative analysis of MR-PET studies with functionalities such as SUV calculation and PET-specific segmentation. In conjunction with the MR Oncology possibilities, multi-modal evaluations and follow-up analyses of MR-PET examinations receive optimum support. |
| 111. | syngo.MR General Engine<br>syngo.MR General Engine extends syngo.via by adding dedicated workflows and tools for routine and advanced reading of MR examinations.                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| 112. | syngo.MR General Engine<br>syngo.MR General Engine extends syngo.via by adding dedicated workflows and tools for routine and advanced reading of MR examinations.                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| 113. | syngo.MR General Engine<br>syngo.MR General Engine extends syngo.via by adding dedicated workflows and tools for routine and advanced reading of MR examinations.                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| 114. | PET Dynamic #mMR<br>Support for offline histogramming and reconstruction. Support for retrospective histogramming in any arbitrary frame durations of 3 second or greater, maximum of 100 frames defined by available disk space.<br>PET Dynamic MRWP #mMR<br>Support for offline histogramming and reconstruction. Support for retrospective histogramming in any arbitrary frame durations of 3 second or greater, maximum of 100 frames defined by available disk space.                                                                                                                                |
| 115. | PET Cardiac #mMR<br>Provides PET cardiac offline histogramming and reconstruction for improved accuracy in quantitation as well as visualization of cardiac motion. Supports a maximum of 16 gate bins from the list mode PET acquisition.<br>PET Cardiac MRWP #mMR<br>Provides PET cardiac offline histogramming and reconstruction for improved accuracy in quantitation as well as visualization of cardiac motion. Supports a maximum of 16 gate bins from the list mode PET acquisition.                                                                                                              |

|             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|-------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>116.</b> | <p>PET Respiratory #mMR<br/>Provides PET respiratory offline histogramming and reconstruction for improved accuracy in quantitation.</p> <p>PET Respiratory MRWP #mMR<br/>Provides PET respiratory offline histogramming and reconstruction for improved accuracy in quantitation.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>117.</b> | <p>Biograph mMR PET Sources, 1st Fill<br/>Calibration sources for the Biograph mMR. These sources are to be purchased with a new Biograph mMR scanner.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>118.</b> | <p>mMR Replacement PET Sources kit<br/>Replacement Calibration sources kit for the Biograph mMR.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>119.</b> | <p>mMR Replacement QC Phantom kit<br/>Replacement QC Phantom only, for the Biograph mMR.<br/>Includes:<br/>Uniform QC Phantom, Ge-68, specific for the Biograph mMR. Max 3.0 mCi (111 MBq)</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>120.</b> | <p>mMR Replacement Line Source kit<br/>Replacement Line Sources only, for the Biograph mMR.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>121.</b> | <p>NEMA PET Self-testing Kit #mMR<br/>The NEMA PET Self-testing Kit includes versions of the following phantoms and fixtures appropriate for use with the Biograph mMR:</p> <ul style="list-style-type: none"> <li>- PET Scatter Phantom NEMA, a 70 cm solid cylinder with a fillable line source</li> <li>- PET Sensitivity Phantom NEMA, a fillable line source and five concentric cylindrical sleeves</li> <li>- PET Image Quality Phantom, a fillable phantom with six fillable spheres and a fillable cylindrical insert</li> <li>- L-bracket for positioning point sources</li> </ul> <p>Two combination holders for positioning NEMA phantoms</p> <ul style="list-style-type: none"> <li>- PET IQ Phantom holder</li> </ul> <p>ACRIN PET Self-testing Kit #mMR<br/>The ACRIN PET Self-testing Kit includes versions of the following phantoms and fixtures appropriate for use with the Biograph mMR:</p> <ul style="list-style-type: none"> <li>- PET ACR Phantom, a fillable phantom with fillable and solid inserts</li> <li>- PET ACR Phantom holder</li> </ul> |

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 122. | Biograph mMR PET Sources, 1st Fill<br>Calibration sources for the Biograph mMR. These sources are to be purchased with a new Biograph mMR scanner.<br>Kit includes:<br>- 4 Ge-68 Line Sources, specific for the Biograph mMR. Max. 1.5 mCi (55.5 MBq) per rod source<br>- Uniform QC Phantom, Ge-68, specific for the Biograph mMR. Max 3.0 mCi (111 MBq)                                                                         |
| 123. | syngo BreVis Biopsy #Tim<br>syngo BreVis Biopsy is a task card for easy and effective breast biopsy planning for the Acquisition Workplace (AWP).<br>syngo BreVis Biopsy MRWP #Tim<br>syngo BreVis Biopsy is a task card for easy and effective breast biopsy planning.<br>syngo BreVis MRWP #Tim<br>syngo BreVis is a task card for optimized visualization of multimodal breast images and for processing dynamic 3D data sets. |
| 124. | Breast Dot Engine #T+D<br>A set of predefined Breast Dot Engine configurations is provided for tissue depiction, implant evaluation and breast biopsy. These configurations offer a comprehensive user guidance, workflow automation and personalization towards the individual patient so that excellent image quality and increased consistency for breast diagnosis and intervention can be achieved.                          |
| 125. | Trigger Converter<br>An optical trigger signal is available to trigger external stimulation devices in fMRI experiments.                                                                                                                                                                                                                                                                                                          |
| 126. | Patient Supervision TV #V,mMR<br>Special video camera for observing the patient during an MR examination, conveniently mounted on the wall.                                                                                                                                                                                                                                                                                       |
| 127. | Patient TV wall support<br>Wall mount for the patient monitor.                                                                                                                                                                                                                                                                                                                                                                    |
| 128. | LCD Color Remote Monitor #T+D<br>High-resolution 19" LCD color flat screen monitor, 1280 x 1024 pixels, with integrated gamma correction for optimum display of MR grayscale images and automatic backlight control for long-term brightness stability.                                                                                                                                                                           |
| 129. | Video Distributor for monitor #T+D<br>Video distributor for connection of up to two additional color LCD monitors as well as an In-Room Acq WP syngo (in the examination room parallel to the syngo Acquisition Workplace operator monitor).                                                                                                                                                                                      |
| 130. | Additional Camera #T+D<br>Additional video camera for monitoring the situation in the examination or in the waiting room (in combination with wide angle option).<br>The additional camera is conveniently mounted on the wall.                                                                                                                                                                                                   |
| 131. | Supervision Camera & Controller #T+D<br>Additional video camera to allow for comprehensive monitoring of the patient's face independent of patient position during an MR examination.                                                                                                                                                                                                                                             |
| 132. | Wide Angle Lens #T+D<br>Wide angle lens to allow for observing the situation in the waiting room.                                                                                                                                                                                                                                                                                                                                 |

|             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|-------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>133.</b> | Application cushion                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>134.</b> | Respiratory cushion with pressure hose                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>135.</b> | Respiratory belt                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>136.</b> | PMU-Display #mMR<br>LCD monitor arm-mounted on the magnet cover for display of physiologic (ECG, pulse, respiration) signals in the examination room.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>137.</b> | Pure White Design #mMR<br>The color of the main face plate cover with integrated control panel and table display is brilliant white surrounded by a brilliant silver trim. The asymmetrical deco area on the left side is colored white matte and also with a brilliant surrounding silver trim.                                                                                                                                                                                                                                                                                                                                                                       |
| <b>138.</b> | Pure White Design #mMR<br>The color of the main face plate cover with integrated control panel and table display is brilliant white surrounded by a brilliant silver trim. The asymmetrical deco area on the left side is colored white matte and also with a brilliant surrounding silver trim.                                                                                                                                                                                                                                                                                                                                                                       |
| <b>139.</b> | Accessories cabinet<br>Specially designed non-ferromagnetic cabinet for easy storage of some of the most commonly used coils and accessories.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>140.</b> | Coil Storage Cart #mMR<br>Specially designed non-ferromagnetic cart for easy storage of the most commonly used coils and accessories.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>141.</b> | Cushions and positioning aids:<br><ul style="list-style-type: none"> <li>- table pad 1:</li> <li>- table pad 2</li> <li>- table pad 3, approved for MR-PET</li> <li>- table pad 4, approved for MR-PET</li> <li>- table pad 5</li> <li>- head-heel insert</li> <li>- paper roll holder, 08112224</li> <li>- 4 belt straps, approved for MR-PET</li> </ul><br><ul style="list-style-type: none"> <li>- table pad 1:<br/>Weight: 1.2 kg<br/>Dimensions: 800 mm x 490 mm x 20 mm</li> <li>- table pad 2<br/>Weight: 1.8 kg<br/>Dimensions: 1100 x 490 x 20 mm</li> <li>- table pad 3, approved for MR-PET<br/>Weight: 0.9 kg<br/>Dimensions: 570 x 490 x 53 mm</li> </ul> |

|             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|-------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|             | <ul style="list-style-type: none"> <li>- table pad 4, approved for MR-PET<br/>Weight: 1.3 kg<br/>Dimensions: 625 x 490 x 53 mm</li> <li>- table pad 5<br/>Weight: 0.3 kg<br/>Dimensions: 490 x 129 x 53 mm</li> <li>- head-heel insert<br/>Weight: 0.945 kg<br/>Dimensions: 300 x 254 x 50 mm</li> <li>- paper roll holder, 08112224<br/>Weight: 1,16 Kg<br/>Dimensions: 540x220x140 mm</li> <li>- 4 belt straps, approved for MR-PET<br/>Weight: 0.1 kg<br/>Dimensions: 1500 x 50 mm</li> </ul> |
| <b>142.</b> | <p>Set of positioning aids</p> <ul style="list-style-type: none"> <li>- knee wedge</li> <li>- positioning roll, approved for MR-PET</li> <li>- back-of-the-head support</li> <li>- 2 large triangles, approved for MR-PET</li> <li>- 2 small triangles, approved for MR-PET</li> <li>- wedge, approved for MR-PET</li> </ul>                                                                                                                                                                     |

|             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|-------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|             | <ul style="list-style-type: none"> <li>- lordosis cushion, approved for MR-PET</li> <li>- knee wedge:<br/>Weight 1.25 kg<br/>Dimensions 557 × 395 × 165 mm</li> <li>- positioning roll, approved for MR-PET:<br/>Weight 1.4 kg<br/>Dimensions 220 × 220 × 399 mm</li> <li>- back-of-the-head support:<br/>Weight 0.45 kg<br/>Dimensions 250 × 250 × 10 mm</li> <li>- 2 large triangles, approved for MR-PET:<br/>Weight 0.4 kg<br/>Dimensions 235 × 176 × 117 mm</li> <li>- 2 small triangles, approved for MR-PET:<br/>Weight 0.4 kg<br/>Dimensions 176 × 94 × 48 mm</li> <li>- wedge, approved for MR-PET:<br/>Weight 0.35 kg<br/>Dimensions 237 × 176 × 67 mm</li> <li>- lordosis cushion, approved for MR-PET:<br/>Weight 0.4 kg<br/>Dimensions 345 × 247 × 40 mm</li> </ul> |
| <b>143.</b> | <p>mMR Head/Neck accessories:</p> <ul style="list-style-type: none"> <li>- Relax cushion</li> <li>- head positioning wedge</li> <li>- cushion neck, approved for MR-PET</li> <li>- head support, short</li> <li>- mirror, approved for MR-PET</li> <li>- head support, approved for MR-PET</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |

|             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|-------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|             | <ul style="list-style-type: none"> <li>- Relax cushion:<br/>Weight 0.4 kg<br/>Dimensions 310 × 250 × 80 mm</li> <li>- head positioning wedge:<br/>Weight 0.05 kg<br/>Dimensions 140 × 80 × 75 mm</li> <li>- cushion neck, approved for MR-PET:<br/>Weight 0.1 kg<br/>Dimensions 290 × 213 × 60 mm</li> <li>- head support, short:<br/>Weight: 0.1 kg<br/>Dimensions: 140 x 100 x 50 mm</li> <li>- mirror, approved for MR-PET:<br/>Weight: 0.05 kg<br/>Dimensions: 132.6 x 100 x 74.3 mm</li> <li>- head support, approved for MR-PET:<br/>Weight: 0.150 kg<br/>Dimensions: 247mm x 247 mm x 98 mm</li> </ul> |
| <b>144.</b> | <p>mMR Spine accessories:</p> <ul style="list-style-type: none"> <li>- cushion Spine matrix, approved for MR-PET<br/>Weight: 2.7 kg<br/>Dimensions: 1185 x 55 x 18 mm</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>145.</b> | <p>mMR Body accessories:</p> <ul style="list-style-type: none"> <li>- mMR Adapter:<br/>Weight: 1.485 kg<br/>Dimensions: 133 x 112.4 x 102.3 mm</li> <li>- 2 straps, approved for MR-PET:<br/>Weight: 0.05 kg<br/>Dimensions: 480 x 50 mm</li> <li>- 4 belts, approved for MR-PET:<br/>Weight: 0.1 kg<br/>Dimensions: 2000 x 50 mm</li> </ul>                                                                                                                                                                                                                                                                  |
| <b>146.</b> | <p>Flex/Loop accessories:</p> <ul style="list-style-type: none"> <li>- hand wedge</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |

|             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|-------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|             | <ul style="list-style-type: none"> <li>- 4 belts</li> <li>- iPAT cushion, small</li> <li>- iPAT cushion, large</li> <li>- 2 soft straps</li> </ul> <p>- hand wedge:<br/>Weight for Flex Coil, large: 1.0 kg<br/>for Flex Coil, small: 0.3 kg<br/>Dimensions for Flex Coil, large:<br/>504 × 212 × 14 mm<br/>for Flex Coil, small:<br/>354 × 162 × 14 mm</p> <p>- 4 belts<br/>Weight: 0,100 kg<br/>Dimensions: 2000 x 50 mm</p> <p>- iPAT cushion, small:<br/>Weight: 0,030 kg<br/>Dimensions: 354 x 162 x 14 mm</p> <p>- iPAT cushion, large:<br/>Weight: 0,050 kg<br/>Dimensions: 504 x 212 x 14 mm</p> <p>2 soft straps:<br/>Weight 0.11 kg<br/>Dimensions 2000 × 50 × 1 mm</p> |
| <b>147.</b> | <p>mMR Head 32 accessories:</p> <ul style="list-style-type: none"> <li>- mirror carrier with mirror</li> <li>- mMR Head 32 Relax cushion</li> <li>- head alignment wedges</li> </ul> <p>mirror carrier with mirror:<br/>Weight 0.25 kg<br/>Dimensions 186 × 110 × 120 mm</p> <p>mMR Head 32 Relax cushion:<br/>Weight 0.3 kg<br/>Dimensions 230 × 275 × 155 mm</p> <p>head alignment wedges:<br/>Weight 0.03 kg<br/>Dimensions 135 × 60 × 50 mm</p>                                                                                                                                                                                                                               |

- 148.** mMR Breast Coil Accessories:
- insertion plate (for medial access);
  - horizontal and vertical slat plates;
  - needle guide (post/pillar method);
  - height adjustable grid (reusable);
  - reusable grid;
  - height adjustable grid (disposable);
  - frame fixture;
  - head rest;
  - abdominal pad;
  - foot support roller.

Please note that the following biopsy accessories are only for training purposes: insertion plate (for medial access), horizontal and vertical slat plates, needle guide (post/pillar method), height adjustable grid (reusable), reusable grid, height adjustable grid (disposable), frame fixture.

Accessories for the examination (needles, needle guide, height adjustable grid, etc.) need to be ordered separately, for example, from NORAS MRI products GmbH. Here are only unsterile components available; sterile components have to be ordered directly at the supplier, for example, from NORAS MRI products GmbH.

The following positioning aids are available to ensure safe and comfortable patient positioning: head rest, abdominal pad, foot support roller.

Insertion plate (for medial access)

Dimensions: 321 x 191 x 37 mm

Weight: 0,192 kg

Horizontal and vertical slat plates

Horizontal:

Dimensions: 170 x 126 x 17 mm

Weight: 0,046 kg

Vertical:

Dimensions: 170 x 104 x 17 mm

Weight: 0,046 kg

Needle guide (post/pillar method)

Dimensions: 10 x 10 x 40 mm

Weight: 0,002 kg

Height adjustable grid (reusable)

Dimensions: 165 x 103 x 15 mm

Weight: 0,078 kg

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>Reusable grid<br/>Dimensions: 272 x 106 x 30 mm<br/>Weight: 0,192 kg</p> <p>Height adjustable grid (disposable)<br/>Dimensions: 164 x 106 x 15 mm<br/>Weight: 0,063 kg</p> <p>Frame fixture<br/>Dimensions: 522 x 89 x 55 mm<br/>Weight: 0,32 kg</p> <p>Head rest<br/>Dimensions: 238 x 196 x 43 mm<br/>Weight: 0,030 kg</p> <p>Abdominal pad<br/>Dimensions: 547 x 453 x 180 mm<br/>Weight: 0,901 kg</p> <p>Foot support roller<br/>Dimensions: 220 x 400 mm<br/>Weight: 0.451 kg</p> |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**149.** Positioning frame for BreastCoil#mMR with accessories:

- a head cushion
- abdominal support cushions
- breast stabilizing cushions (4 pcs.)

Open design positioning frame for inserting the mMR Breast Coil.

The accessories are available to ensure safe and comfortable patient positioning.

**mMR Frame Breast**

Dimensions: 736 x 483 x 174 mm

Weight: 5,108 kg

**Head Rest Breast**

Dimensions: 238 x 196 x 43 mm

Weight: 0,030 kg

**Cushion Breast**

Dimensions: 547 x 453 x 180 mm

Weight: 0,901 kg

**Breast stabilizing cushions**

**Small:**

Dimensions: 141x 118 x 50 mm

Weight: 0,019 kg

**Large:**

Dimensions: 141 x 118 x 117 mm

Weight: 0,044 kg

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 150. | <p>Positioning frame for BreastCoil with accessories:</p> <ul style="list-style-type: none"> <li>- a head cushion</li> <li>- abdominal support cushions</li> <li>- breast stabilizing cushions (4 pcs.)</li> </ul> <p>Open design positioning frame for inserting the 4-channel BI Breast Coil or 16-Channel AI Breast Coil. The accessories are available to ensure safe and comfortable patient positioning.</p> <p><b>Frame Breast</b><br/>Dimensions: 736 x 483 x 174 mm<br/>Weight: 5,108 kg</p> <p><b>Head Rest Breast</b><br/>Dimensions: 238 x 196 x 43 mm<br/>Weight: 0,030 kg</p> <p><b>Cushion Breast</b><br/>Dimensions: 547 x 453 x 180 mm<br/>Weight: 0,901 kg</p> <p><b>Breast stabilizing cushions</b><br/>Small:<br/>Dimensions: 141x 118 x 50 mm<br/>Weight: 0,019 kg<br/>Large:<br/>Dimensions: 141 x 118 x 117 mm<br/>Weight: 0,044 kg</p> |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**151.** TxRx 15 Ch Knee/TxRx CP Extremity accessories:

- coil pad, 0.25
- coil pad, 0.50
- ankle pad
- other knee pad

The accessories are available to ensure safe and comfortable patient positioning.

Coil pad, 0.25

Dimensions: 138 x 50 x 254 mm

Weight: 0.350 kg

Coil pad, 0.50

Dimensions: 138 x 50 x 254 mm

Weight: 0.224 kg

Ankle pad

Dimensions: 150 x 150 x 300 mm

Weight: 0.6 kg

Other Knee pad

Dimensions: 200 x 172 x 400 mm

Weight: 0.85 kg

**152.** Shoulder Array accessories:

- Cup cushion
- Special head support
- Pad

The accessories are available to ensure safe and comfortable patient positioning.

Cup cushion

Dimensions: 225 x 225 x 13 mm

Weight: 0.35 kg

Special head support

Dimensions: 294 x 225 x 93 mm

Weight: 0.55 kg

Pad

Dimensions: 425 x 345 x 20mm

Weight: 0.9 kg

east Coil.

|             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|-------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>153.</b> | Biopsy Starter Kit #Tim<br>Please note that the scope of delivery from Siemens provides biopsy accessories only for training purposes. Accessories for the examination (needles, needle guide, height adjustable grid, etc.) need to be ordered separately, for example, from NORAS MRI products GmbH. Here are only unsterile components available; sterile components have to be ordered directly at the supplier, for example from NORAS MRI products GmbH. |
| <b>154.</b> | Installation material                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>155.</b> | Footswitch<br>MR foot pedal with two pneumatic buttons for switching on and off the preset MR sequence. The foot pedal is MR compatible and is located next to the patient table in the examination room. The foot pedal is especially convenient for sequence control during MR examinations during surgical procedures                                                                                                                                       |
| <b>156.</b> | PC Keyboard US English #Tim<br>Standard PC keyboard with 101 keys                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>157.</b> | Biograph mMR – System<br>License key for integrated software                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>158.</b> | Video Distributor for monitor #T+D<br>Video distributor for connection of up to two additional color LCD monitors as well as an In-Room Acq WP syngo (in the examination room parallel to the syngo Acquisition Workplace operator monitor).                                                                                                                                                                                                                   |
| <b>159.</b> | UPS system #mMR<br>UPS system EATON 9PX3000 VA for Biograph mMR for safeguarding computers. Including Power Cable of 9 m for connecting the UPS.                                                                                                                                                                                                                                                                                                               |
| <b>160.</b> | UPS Battery module #mMR<br>UPS battery module EATON Eaton 9PX 3kVA for Biograph mMR for safeguarding computers                                                                                                                                                                                                                                                                                                                                                 |
| <b>161.</b> | Distributor MR 160A                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>162.</b> | Distribution board                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>163.</b> | RF-Cabin, Steel #mMR<br>Steel full-room RF cabin, no radiation protection, module method of building, for Biograph mMR. Maximum outside dimensions:<br>Length x Width up to 35 sqm<br>Height up to 3,00 m                                                                                                                                                                                                                                                      |
| <b>164.</b> | Filter plate for RF-Cabin<br>Filterplate in low leakage circuitry for shielded rooms (MRI)<br>6 separate lines, 250 V, 50/60 Hz, 250Vdc, 32A<br>For usage at MR-Shielding system the filterplate is manufactured completely in non-magnetic materials.                                                                                                                                                                                                         |
| <b>165.</b> | RF-Traverse Beams, Steel #3T<br>5 traverse beams for the RF cabin (Art. No. 14413606 and 14441589) of the MAGNETOM Trio a Tim System, Verio, Prisma as well as Biograph mMR.<br>Length 4,590 mm                                                                                                                                                                                                                                                                |

|             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|-------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>166.</b> | <p>RF-Interior #3T</p> <p>Standard interior décor with Halogen lighting for the RF cabin (Art. No. 14413606 and 14441589) of the MAGNETOM Trio a Tim System, Verio, Prisma as well as Biograph mMR.</p> <p>Standard outside dimensions:<br/>L x W x H = 7,590 x 4,590 x 3,000 mm</p> <p>RF-Interior, LED #3T</p> <p>Standard interior décor with LED lighting for the RF cabin (Art. No. 14413606 and 14441589) of the MAGNETOM Trio a Tim System, Verio, Prisma as well as Biograph mMR.</p> <p>Maximum outside dimensions:<br/>Length x Width up to 35 sqm<br/>Height up to 3,00 m</p>                                       |
| <b>167.</b> | <p>RF Cabin Sound Proofing, Steel #3T</p> <p>Sound protection system for MAGNETOM Trio a Tim System, Verio and Prisma steel modular cabin (Art. No. 14413606) as well as Biograph mMR steel modular cabin (Art. Nr. 14441589).</p> <p>Maximum outside dimensions:<br/>Length x Width up to 35 sqm<br/>Height up to 3,00 m</p> <p>Due to building requirements according to DIN 4109 (sound protection for structural engineering) it is possible to get the following values:<br/>Sound protection, walls &gt; 52 dB(A)<br/>Sound protection, ceiling &gt; 56 dB(A)<br/>Sound protection, observation window &gt; 50 dB(A)</p> |
| <b>168.</b> | <p>Filter Cabinet, big</p> <p>Cabinet for filter plate inside of the RF cabin for Biograph mMR. Wall-mounted cabinet custom-made. Outside dimensions appr. 1.300 x 420 x 2.700 mm (L x W x H)</p>                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>169.</b> | <p>Threshold-free door</p> <p>Soundproof threshold-free door for the RF-Cabins of MAGNETOM 1.5T/3T that can be opened and closed without effort by pneumatic mechanism.</p> <p>Clear opening: W x H 1,200 x 2,050 mm</p>                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>170.</b> | <p>Emergency hatch</p> <p>With decorative plate on both sides, white, to be opened from the outside, approx. 600 x 600 mm.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>171.</b> | <p>Observation/ Day-light window</p> <p>Observation / daylight window that can be opened, installation sideways.</p> <p>Dimensions: 1,200 x 800 mm.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>172.</b> | <p>MR Workplace Table 1.2m</p> <p>Table suited for syngo Acquisition Workplace and syngo MR Workplace based on syngo Hardware.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>173.</b> | <p>MR Workplace Container, 50cm</p> <p>50 cm wide extra case for the syngo host computer with sliding front door to allow change of storage media (CD/DVD/USB).</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>174.</b> | <p>Swivel chair</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |

|             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>175.</b> | IV-Rod, non-magnetic f. Stretcher<br>IV-Rod, fix mounted, for AGA nonmagnetic stretcher.<br>- Diameter: 180mm<br>- Length:1720mm<br>- With four hooks                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>176.</b> | Instrument table, non-magnetic<br>Mobile 4-leg frame construction made of 18/10 CrNi steel.<br>Optimal - referring to the handlebar - working height of 800 mm. 63.5 cm x 43.5 cm tubular frame forms a safety rail; no sharp edges. Two shelves: surface ground to reduce glare, resistant to disinfectants, removable. Standard sizes, dimensioned to conform with the DIN standard for sterile goods. 4 plastic swivel castors, two of them with brakes, wear resistant, non-discoloring when in contact with other material, electrically conductive, Ø 75 mm, 4 wall spacer rings. |
| <b>177.</b> | Metal detector                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>178.</b> | Fire extinguisher, anti-magn., emp.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>179.</b> | Emergency switches<br>The MR system has different types of emergency switches.<br>Magnet Stop Shutting down the static magnetic field (quenching) - in case of accidents with attracted metal parts or in case of fire emergency shutdown electric power of the entire MR system is switched off, but magnet remains at field e.g. in case of fire Table Stop Motorized table movement is stopped<br>In case of accidents or injury due to table movement                                                                                                                               |
| <b>180.</b> | Earthquake Kit for ECO Chiller 60 kW.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>181.</b> | ECG-electrodes                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>182.</b> | CD-disk for storage images                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>183.</b> | DVD offer higher storage capacity than compact discs while having the same dimensions.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>184.</b> | Service Syphon Kit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>185.</b> | HE-SYPHON DEWAR END                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>186.</b> | Integrated magnet cooling system                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>187.</b> | Floor plate for table                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>188.</b> | Large Joint Dot Engine E11 #T+D<br>Large Joint Dot Engine optimizes image quality of knee, hip and shoulder scans by proposing the most appropriate protocols according to the examination strategy chosen for the specific patient. It ensures reproducible image quality and streamlines large joint examinations to the greatest extent.                                                                                                                                                                                                                                             |
| <b>189.</b> | Large Joint Dot Engine E11 #T+D<br>Large Joint Dot Engine optimizes image quality of knee, hip and shoulder scans by proposing the most appropriate protocols according to the examination strategy chosen for the specific patient. It ensures reproducible image quality and streamlines large joint examinations to the greatest extent.                                                                                                                                                                                                                                             |

|             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|-------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>190.</b> | <b>Large Joint Dot Engine E11 #T+D</b><br>Large Joint Dot Engine optimizes image quality of knee, hip and shoulder scans by proposing the most appropriate protocols according to the examination strategy chosen for the specific patient. It ensures reproducible image quality and streamlines large joint examinations to the greatest extent.                                                                                                                   |
| <b>191.</b> | <b>Spine Dot Engine #mMR</b><br>The Spine Dot Engine provides optimized cervical, thoracic and lumbar spine imaging. Amongst various features to support a streamlined spine workflow is AutoLabeling of the vertebrae.                                                                                                                                                                                                                                              |
| <b>192.</b> | <b>Advanced Diffusion #mMR</b><br>QuietX DWI and RESOLVE together make up the Advanced Diffusion package.<br>QuietX DWI enables quieter diffusion-weighted imaging of the brain with up to 70% reduction in sound pressure relative to conventional diffusion-weighted imaging.<br>RESOLVE (Readout Segmentation Of Long Variable Echo-trains) is a multi-shot, readout segmented EPI sequence for high-resolution, low-distortion diffusion-weighted imaging (DWI). |
| <b>193.</b> | <b>Breast Dot Engine #T+D</b><br>A set of predefined Breast Dot Engine configurations is provided for tissue depiction, implant evaluation and breast biopsy. These configurations offer a comprehensive user guidance, workflow automation and personalization towards the individual patient so that excellent image quality and increased consistency for breast diagnosis and intervention can be achieved.                                                      |
| <b>194.</b> | <b>Adhes. Plastic-Mastic complete</b>                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>195.</b> | <b>QISS #T+D</b><br>Software package with QISS sequence, protocols and Dot AddIn for non-contrast enhanced peripheral MRA. QISS particularly enables higher reproducibility than existing methods and is an alternative to MR angiography techniques with contrast medium, especially for patients with severe renal insufficiency.                                                                                                                                  |
| <b>196.</b> | <b>FREEZEit Body MRI Package #T+D</b><br>FREEZEit Body Package contains two robust sequences for advanced body imaging: TWIST-VIBE and StarVIBE.                                                                                                                                                                                                                                                                                                                     |
| <b>197.</b> | <b>Breast Biopsy Software #V,mMR,Sk</b><br>This syngo based postprocessing software allows the user to perform the localization and coordinate planning for inserting the needle using the 'post and pillar' method. The user-friendly software computes the coordinates after the center of the lesion and the marker of the breast biopsy unit have been selected.                                                                                                 |
| <b>198.</b> | <b>Seismic Kit #T+D</b>                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>199.</b> | <b>PMU-Display #mMR</b><br>LCD monitor arm-mounted on the magnet cover for display of physiologic (ECG, pulse, respiration) signals in the examination room.                                                                                                                                                                                                                                                                                                         |
| <b>200.</b> | <b>Pure White Design #mMR</b><br>The color of the main face plate cover with integrated control panel and table display is brilliant white surrounded by a brilliant silver trim. The asymmetrical deco area on the left side is colored white matte and also with a brilliant surrounding silver trim.                                                                                                                                                              |

|             |                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|-------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>201.</b> | <p>Body Compass #mMR</p> <p>BodyCompass software enables periodical respiratory motions of the patient to be acquired and modeled in the volume to be examined by means of MR imaging, using a radial VIBE sequence.</p>                                                                                                                                                                                                                    |
| <b>202.</b> | <p>Advanced WARP #T+D</p> <p>Advanced WARP application consists of SEMAC, a technique to reduce gross metal artifacts (i.e. through-plane artifacts) caused by big orthopedic implants. The main clinical applications are in hip and knee joint replacements. Available protocols include T1-weighted, T2-weighted, proton density and STIR contrast.</p>                                                                                  |
| <b>203.</b> | <p>Integrated magnet cooling system</p>                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>204.</b> | <p>Eco Chiller 60kW</p> <p>The KKT ECO 133 - L chiller is a dedicated 20°C cooling system which automatically adapts to the different cooling requirements (e.g. system in operation, standby, ...) to reduce the energy consumption for cooling.</p> <p>The cooling system must be used in combination with the IFP (Interface Panel), if there is no on-site chilled water supply at all. The IFP is included in the scope of supply.</p> |

### 5.3 Technical Requirements

Technical Requirements based on the key parameters and values as presented in the chapter below.

#### ■ Gradient data

|                                                                      |             |
|----------------------------------------------------------------------|-------------|
| MQ Gradients                                                         |             |
| Performance for each axis                                            |             |
| Max. amplitude <sup>1)</sup>                                         | 45 mT/m     |
| Min. rise time                                                       | 225 $\mu$ s |
| Max. slew rate                                                       | 200 T/m/s   |
| Vector gradient performance (vector addition of all 3 gradient axes) |             |
| Max. eff. amplitude                                                  | 78 mT/m     |
| Max. eff. slew rate                                                  | 346 T/m/s   |

**1) Values is valid for each direction**

#### ■ Magnet system

Short-bore, patient-friendly design, high homogeneity

Active Shielding (AS) for easy installation

Integrated External Interference (E.I.S.) shielding provides permanent protection against potential interferences through ferromagnetic objects moved into the proximity of the magnet

Operating mode: First Level controlled operating mode according to IEC 60601-2-33.

| Magnet data                             |                     |
|-----------------------------------------|---------------------|
| Operating field strength                | 3 Tesla             |
| Magnet type                             | Niobium-Titanium    |
| Field stability over time               | <0.1 ppm/h          |
| Weight (with cryogenics)                | 6300 kg             |
| Magnet length                           | 163 cm              |
| Bore size                               | 60 cm               |
| System length cover to cover            | 199 cm              |
| Magnet Cooling system                   |                     |
| Refill interval (typical) <sup>1)</sup> | 10 years            |
| Boil-off rate (typical) <sup>1)</sup>   | 0.0 l/year          |
| Max. helium capacity                    | approx. 1500 liters |
| Minimum helium level                    | 50%                 |
| Cryostat                                | Stainless steel     |

**1) For typical clinical use, depending on sequences and operating time with running helium compressor.**

**The system needs to be serviced at regular interval. Undisturbed magnet cooling for 24 hours and 7 days a week.**

| Magnetic field homogeneity   |           |
|------------------------------|-----------|
| 45 x 50 x 50 cm <sup>3</sup> | ≤4.0 ppm  |
| 40 cm DSV                    | ≤1.4 ppm  |
| 30 cm DSV                    | ≤0.3 ppm  |
| 20 cm DSV                    | ≤0.05 ppm |
| 10 cm DSV                    | ≤0.01 ppm |

DSV = Diameter spherical volume (x, y and z direction)

Standard deviation  $V_{rms}$  (Volume root-mean square) measured with highly accurate 24 plane plot method (20 per plane).

Standard active shim with 3 linear and 5 non-linear channels (2nd order)

| Magnet field shimming                                                                                                                             |              |
|---------------------------------------------------------------------------------------------------------------------------------------------------|--------------|
| <ul style="list-style-type: none"> <li>• Patient-specific automated shim</li> <li>• For optimization of the magnetic field homogeneity</li> </ul> |              |
| Time to shim                                                                                                                                      | approx. 20 s |

| Magnet shielding            |                                                                    |             |
|-----------------------------|--------------------------------------------------------------------|-------------|
| Active Shielding (AS)       | 5th generation active shielding (AS) technology with counter coils |             |
| Fringe field (axial×radial) | 0.5 mT <sup>1)</sup>                                               | 4.6 m×2.6 m |
|                             | 0.1 mT                                                             | 6.8 m×3.9 m |

1) pacemaker safety limit

#### ■ Patient handling

| MR Tim Table                                      |                      |                  |
|---------------------------------------------------|----------------------|------------------|
| Max. patient weight (including vertical movement) |                      | 200 kg           |
| Max. scan range                                   |                      | 200 cm           |
| Vertical table movement                           | Range                | approx. 50–90 cm |
|                                                   | Max. speed           | 10 cm/s          |
| Horizontal table movement                         | Max. range           | 291.0 cm         |
|                                                   | Max. speed           | 20 cm/s          |
|                                                   | Positioning accuracy | ±0.8 mm          |

#### ■ Hearing protection data

The A-evaluated, effective sound pressure level was measured according to NEMA MS 4-2006 (National Electrical Manufacturers Association) using the maximum gradient acoustic noise (MGAN) method.

#### Patient noise

Patients require hearing protection with an  $SNR^1=23$  dB or more.

<sup>1</sup> SNR = Single Number Rating

**Noise for personnel in the examination room**

Noise measured: 87.9 dB (A)

Hearing protection: SNR=6 dB or more

▪ **MR gradient system (MQ gradients)**

| General features                                                                                               |            |
|----------------------------------------------------------------------------------------------------------------|------------|
| Actively shielded (AS) whole-body gradient coil system                                                         |            |
| Extremely low eddy currents                                                                                    |            |
| Water-cooled coil and amplifier for maximum performance                                                        |            |
| Gradient coil length 159 cm, optimized for gradient linearity                                                  |            |
| Excellent spatial and temporal resolution                                                                      |            |
| Shortest TE and TR for high image contrast and best signal-to-noise ratio                                      |            |
| Patient-specific shimming                                                                                      |            |
| Active shim via 3 linear channels with 20 coils                                                                |            |
| Combined with 5 additional non-linear 2nd order channels with 32 coils for advanced high order shim (standard) |            |
| Maximum acoustic noise level                                                                                   | 118 dB (A) |

| Gradient amplifier                                               |        |
|------------------------------------------------------------------|--------|
| Water-cooled, highly compact, modular design                     |        |
| Ultra-fast solid-state technology with very low switching losses |        |
| Max. output voltage <sup>1)</sup>                                | 2250 V |
| Max. output current <sup>1)</sup>                                | 900 A  |

**1) Values for each of the 3 gradient axes**

| Gradient performance                                                 |           |
|----------------------------------------------------------------------|-----------|
| Performance for each axis                                            |           |
| Max. amplitude <sup>1)</sup>                                         | 45 mT/m   |
| Min. rise time                                                       | 225 µs    |
| Max. slew rate                                                       | 200 T/m/s |
| Vector gradient performance (vector addition of all 3 gradient axes) |           |
| Max. eff. amplitude                                                  | 78 mT/m   |
| Max. eff. slew rate                                                  | 346 T/m/s |

**1) Values is valid for each direction**

| Resolution parameters |      |                      |
|-----------------------|------|----------------------|
|                       | min  | max                  |
| FOV                   | 5 mm | 500 mm <sup>1)</sup> |

| Resolution parameters       |         |           |
|-----------------------------|---------|-----------|
| Slice thickness 2D          | 0.1 mm  | 200 mm    |
| Partition thickness 3D      | 0.05 mm | 20 mm     |
| Slab thickness 3D           | 5 mm    | 500 mm    |
| Max. matrix                 |         | 1024      |
| Highest in-plane resolution |         | 9 $\mu$ m |

1) Depending on the application and sequence, the maximum FOV in the z-direction can be up to 45 cm

#### RF technology

| RF data                             |                                                              |              |
|-------------------------------------|--------------------------------------------------------------|--------------|
| Peak power of transmitter amplifier | 35 kW                                                        |              |
| Transmitter bandwidth               | 500 kHz                                                      |              |
| Receiver bandwidth                  | 500 Hz – 1 MHz (for each channel)                            |              |
| RF transmit coils                   |                                                              |              |
| Body coil                           | Max. applied RF field B1+                                    | 23.5 $\mu$ T |
|                                     | Max. specified B1+ rms                                       | 3.2 $\mu$ T  |
|                                     | Distance to isocenter, where RF transmit field is reduced by |              |
|                                     | 3 dB                                                         | 0.15 m       |
|                                     | 10 dB                                                        | 0.27 m       |

#### MR digital signal processing system

| General                                              |                         |
|------------------------------------------------------|-------------------------|
| Compact air cooled amplifier                         |                         |
| Digital quadrature demodulation                      |                         |
| Digital filtering                                    |                         |
| Receiver with high dynamic range without adjustments |                         |
| Frequency stability (5 min)                          | $\pm 2 \times 10^{-10}$ |
| Frequency control                                    | 32 bits (0.03 Hz)       |
| Phase control                                        | 16 bits (0.006 degrees) |

| Transmitter path       |                                                                     |
|------------------------|---------------------------------------------------------------------|
| Transmit amplitude     | 16 bit control<br>50 ns resolution                                  |
| Gain stability (5 min) | $\leq 0.1$ dB (low signal path)<br>$\leq 0.15$ dB (power amplifier) |

| Transmitter amplifier        |              |
|------------------------------|--------------|
| Transmit amplifier bandwidth | 800 kHz      |
| Peak power                   | 35 kW        |
| Max. applied RF-Field        | 23.5 $\mu$ T |

| Receive path                                          |                                                                                                |
|-------------------------------------------------------|------------------------------------------------------------------------------------------------|
| Receiver bandwidth                                    | 500 Hz–1 MHz (for each channel)                                                                |
| Receiver signal resolution                            | 32 bit                                                                                         |
| Number of independent receiver channels               | 18/32                                                                                          |
| Max. number of simultaneous quadrature receiver       | 32                                                                                             |
| Sampling rate (sampling resolution)                   | 10 MHz (100 ns)                                                                                |
| Preamplifier noise figure                             | $\leq 0.5$ dB                                                                                  |
| Preamplifier bandwidth                                | 665 kHz                                                                                        |
| Total gain                                            | 31 dB/46 dB with automatic control                                                             |
| Dynamic range (referred to 1 Hz resolution bandwidth) | 162 dB instantaneous at receiver<br>169 dB with automatic gain control at local coil connector |

#### Multi-Nuclear Option<sup>1)</sup>

Multi-Nuclear Imaging and Spectroscopy for the following nuclei:

3He, 7Li, 13C, 17O, 19F, 23Na, 31P, 129Xe

Decoupling nuclei: 13C, 23Na, 31P

Hyperpolarization nuclei: 3He, 13C, 129Xe

| Frequency band                     | 3He      | 7Li      | 13C      | 17O      | 19F      | 23Na     | 31P      | 129Xe    |
|------------------------------------|----------|----------|----------|----------|----------|----------|----------|----------|
| Center frequency [MHz]             | 93.8524  | 47.8799  | 30.9806  | 16.7012  | 115.9040 | 32.5885  | 49.8719  | 34.0816  |
| Usable transmitter bandwidth [kHz] | $\pm 50$ | $\pm 50$ | $\pm 50$ | $\pm 50$ | $\pm 50$ | $\pm 50$ | $\pm 50$ | $\pm 50$ |

| Multi-Nuclear Option <sup>1)</sup> |      |      |     |     |      |     |      |  |
|------------------------------------|------|------|-----|-----|------|-----|------|--|
| Usable receiver bandwidth [kHz]    | ±190 | ±100 | ±65 | ±50 | ±235 | ±65 | ±100 |  |

1) Optional

▪ **PET detector details**

| Detector assembly                |               |
|----------------------------------|---------------|
| Crystal material                 | LSO           |
| Crystal element dimension        | 4 x 4 x 20 mm |
| Crystal elements per block       | 64            |
| Avalanche Photodiodes (APDs)     | 9 per block   |
| Detector ring diameter           | 656 mm        |
| Blocks per detector ring         | 56            |
| Transaxial FOV                   | 588 mm        |
| Plane spacing                    | 2.0 mm        |
| Crystal elements per ring        | 448           |
| Number of detector element rings | 64            |
| Total detector blocks            | 448           |
| Total number of crystal elements | 28672         |
| Axial FOV                        | 258 mm        |
| Image planes                     | 127           |

| Data acquisition/processing |                                                                                                                                              |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|
| System energy resolution    | ≤14.5% FWHM                                                                                                                                  |
| Acquisition modes           | Sinogram, Listmode                                                                                                                           |
| PET Acquisition System      | Supports specialized hardware and software required by the system to acquire and process PET data                                            |
| Standard software           | 3D acquisition<br>Sinogram (static beds)<br>Automated multi-bed acquisitions<br>Standalone PET acquisitions can be performed in service mode |
| Optional software           | Listmode acquisition with offline histogramming supports dynamic and gated im                                                                |

| Performance                                                                                                                                              |
|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| Performance measures represent values derived from internal testing (typical).<br>All measurements are performed with the factory LLD setting of 430 keV |

| Performance |                                                                              |                   |     |
|-------------|------------------------------------------------------------------------------|-------------------|-----|
| NEMA 2012   | Spatial resolution, standard processing (344x344) – axial                    | FWHM @ 1 cm (mm)  | 4.5 |
|             |                                                                              | FWHM @ 10 cm (mm) | 6.8 |
|             | Spatial resolution – transverse (average of radial and tangential) (344x344) | FWHM @ 1 cm (mm)  | 4.2 |
|             |                                                                              | FWHM @ 10 cm (mm) | 5.1 |
|             | Scatter fraction (at peak NEC rate)                                          | 38%               |     |
|             | Sensitivity (cps/kBq)                                                        | 14.2              |     |
|             | Peak NEC rate (kcps)<br>Measured $\leq 26$ kBq/cc                            | 180               |     |

| Additional PET features             |                                                                                                                                  |
|-------------------------------------|----------------------------------------------------------------------------------------------------------------------------------|
| PHS Accessories (PET)<br>(Standard) | Uniform phantom for normalization<br>PET phantom holders<br>Cross-calibration water phantom<br>Alignment phantom + 4 rod sources |

▪ **Data in accordance with IEC 61675-1**

|                                                                                                            |                                                                                     |
|------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| COINCIDENCE WINDOW                                                                                         | 5.56 ns                                                                             |
| Axial sampling                                                                                             | 2.0 mm                                                                              |
| Septal length                                                                                              | Not applicable                                                                      |
| Septal thickness                                                                                           | Not applicable                                                                      |
| Type of transmission source and source ACTIVITY (nominal and recommended range)                            | Not applicable                                                                      |
| Any additional information being considered essential by the MANUFACTURER to characterize normal operation | Not applicable                                                                      |
| AXIAL SLICE WIDTH according to 3.1.5                                                                       | 4.6 mm FWHM                                                                         |
| axial step size according to 3.1.5                                                                         | 2.0 mm                                                                              |
| NORMALIZED SLICE SENSITIVITY                                                                               | 14.1 Sensitivity (cps/kBq) / 127 slices<br>0.11 cps/kBq (Average slice sensitivity) |
| measured TRUE COINCIDENCE COUNT RATE at 20 % COUNT LOSS and 50 % COUNT LOSS                                | at 20%: $563 \times 10^3$ cps<br>at 50%: $617 \times 10^3$ cps                      |
| ADDRESS PILE UP                                                                                            | Not significant                                                                     |

▪ **Coils**

The following standard coils are optimized for minimized PET attenuation and can be used for combined MR and PET exams.

| Coil name                                           | Applications                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |         |        |                         |                                          |
|-----------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------|--------|-------------------------|------------------------------------------|
| mMR Head/Neck, A Tim coil<br>(basic and additional) | <ul style="list-style-type: none"> <li>• Head or neck examinations</li> <li>• Combined head / neck examination</li> <li>• TMJ (temporo mandibular joints)</li> </ul> <p>16-element coil design with 16 integrated preamplifiers, one ring of 12 elements (4 clusters of 3 elements each) and one ring of 4 elements (2 clusters of 2 elements each).</p> <ul style="list-style-type: none"> <li>• 4-channel CP Mode</li> <li>• 8-channel Dual Mode</li> <li>• 16-channel Triple Mode</li> </ul> <p>Pure receiver coil (Rx)</p> <table> <tr> <td>Weight:</td><td>5 kg</td></tr> <tr> <td>Dimensions (L × W × H):</td><td>470 × 300 × 330 mm</td></tr> </table> | Weight: | 5 kg   | Dimensions (L × W × H): | 470 × 300 × 330 mm                       |
| Weight:                                             | 5 kg                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |         |        |                         |                                          |
| Dimensions (L × W × H):                             | 470 × 300 × 330 mm                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |         |        |                         |                                          |
| mMR Head 32, A 3T Tim Coil                          | <ul style="list-style-type: none"> <li>• Head examination</li> </ul> <p>32-element coil design with 32 integrated preamplifiers.<br/>Attenuation correction for simultaneous PET-MR imaging.<br/>For PET-MR usage, both coil parts need to be mounted.</p> <p>Pure receiver coil (Rx)</p> <table> <tr> <td>Weight:</td><td>5.6 kg</td></tr> <tr> <td>Dimensions (L × W × H):</td><td>330 × 374 × 298 mm (Coil without cables)</td></tr> </table>                                                                                                                                                                                                              | Weight: | 5.6 kg | Dimensions (L × W × H): | 330 × 374 × 298 mm (Coil without cables) |
| Weight:                                             | 5.6 kg                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |         |        |                         |                                          |
| Dimensions (L × W × H):                             | 330 × 374 × 298 mm (Coil without cables)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |         |        |                         |                                          |
| mMR Spine, A Tim coil (basic and additional)        | <ul style="list-style-type: none"> <li>• High resolution imaging of the whole spine</li> </ul> <p>Various applications in combination with additional coils</p> <p>24-element coil design with 24 integrated preamplifiers, 8 clusters of 3 elements each.</p> <ul style="list-style-type: none"> <li>• 8-channel CP Mode</li> <li>• 16-channel Dual Mode</li> <li>• 24-channel Triple Mode</li> </ul> <p>Pure receiver coil (Rx)</p> <table> <tr> <td>Weight:</td><td>11 kg</td></tr> <tr> <td>Dimensions (L × W × H):</td><td>1185 × 485 × 33 mm</td></tr> </table>                                                                                         | Weight: | 11 kg  | Dimensions (L × W × H): | 1185 × 485 × 33 mm                       |
| Weight:                                             | 11 kg                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |         |        |                         |                                          |
| Dimensions (L × W × H):                             | 1185 × 485 × 33 mm                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |         |        |                         |                                          |

| Coil name                                         | Applications                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |                   |
|---------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------|
| mMR Body, A Tim coil (basic and additional)       | <ul style="list-style-type: none"> <li>• Thorax (incl. heart)</li> <li>• Abdomen</li> <li>• Pelvis</li> <li>• Hip</li> <li>• Whole-body</li> </ul> <p>6-element coil design with 6 integrated preamplifiers, with 2 clusters of 3 elements each.</p> <ul style="list-style-type: none"> <li>• 2-channel CP Mode</li> <li>• 4-channel Dual Mode</li> <li>• 6-channel Triple Mode</li> </ul> <p>which results in 4-, 8-, or 12-channel body imaging in combination with the Spine Matrix.</p> <p>Pure receiver coil (Rx)</p> |                   |
|                                                   | Weight:                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    | 1.5 kg            |
|                                                   | Dimensions (L × W × H):                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    | 322 × 520 × 70 mm |
|                                                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |                   |
| 4-Channel Flex coil, large (basic and additional) | <ul style="list-style-type: none"> <li>• Imaging of large regions such as medium to large shoulder, hip, and knee</li> </ul> <p>Four integrated low-noise preamplifiers.</p> <p>Pure receiver coil (Rx)</p>                                                                                                                                                                                                                                                                                                                |                   |
|                                                   | Weight:                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    | 0.55 kg           |
|                                                   | Dimensions (L × W × H):                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    | 516 × 224 mm      |
| 4-Channel Flex coil, small (basic and additional) | <p>Imaging of small regions such as small to medium shoulders, wrist, elbow and ankle</p> <p>Four integrated low-noise preamplifiers.</p> <p>Pure receiver coil (Rx)</p>                                                                                                                                                                                                                                                                                                                                                   |                   |
|                                                   | Weight:                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    | 0.45 kg           |
|                                                   | Dimensions (L × W × H):                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    | 366 × 174 mm      |

|                           |                                                                                                                                                                                                                                                                                                                                          |                    |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------|
| mMR Breast Coil           | <p>The coil consists of a 2-channel insertion plate and a 2-channel support coil each and is optimized for the combined MR/PET mode.</p> <p>Breast imaging.<br/>4 Channel phased array.<br/>Receive.<br/>Simultaneous basic imaging of both breasts in all directions.<br/>Uni-lateral biopsy imaging for lateral and medial access.</p> |                    |
|                           | Weight:                                                                                                                                                                                                                                                                                                                                  | 5.5 kg             |
|                           | Dimensions (L x W x H):                                                                                                                                                                                                                                                                                                                  | 525 x 379 x 40 mm  |
| Tx/Rx CP Extremity Coil   | <p>Transmit/receive coil for joint examinations in the region of the lower extremities.</p> <p>Examinations of joints in the area of the lower extremities, e.g. ankle and knee.<br/>1-Channel Circular-Polarized (CP) Volume Coil<br/>Transmit/Receive.</p>                                                                             |                    |
|                           | Weight:                                                                                                                                                                                                                                                                                                                                  | 8 kg               |
|                           | Dimensions (L x W x H):                                                                                                                                                                                                                                                                                                                  | 532 x 315 x 419 mm |
| Small Shoulder Array Coil | <p>Receive array coil, anatomically adapted coil tops attached to a base plate, either for the left or for the right shoulder.</p> <p>Examination of the shoulder.<br/>4-Channel coil design with 4 integrated preamplifiers.<br/>Receive.</p>                                                                                           |                    |
|                           | Weight:                                                                                                                                                                                                                                                                                                                                  | 3 kg               |
|                           | Dimensions (L x W x H):                                                                                                                                                                                                                                                                                                                  | 185 x 200 x 190 mm |
|                           | <p>Base plate</p> <p>For optimal image quality with slim patients, the shoulder coil cups should be positioned at the inside guide rails of the base plate. The outer guide rails should be used only for muscular and adipose patients.</p>                                                                                             |                    |
|                           | Weight:                                                                                                                                                                                                                                                                                                                                  | 5.4 kg             |
|                           | Dimensions (L x W x H):                                                                                                                                                                                                                                                                                                                  | 489 x 444 x 46 mm  |

|                           |                                                                                                                                                                                                                     |                    |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------|
| Large Shoulder Array Coil | Receive array coil, anatomically adapted coil tops attached to a base plate, either for the left or for the right shoulder.                                                                                         |                    |
|                           | Examination of the shoulder.<br>4-Channel coil design with 4 integrated preamplifiers.<br>Receive.                                                                                                                  |                    |
|                           | Weight:                                                                                                                                                                                                             | 3 kg               |
|                           | Dimensions (L × W × H):                                                                                                                                                                                             | 220 x 210 x 190 mm |
|                           | Base plate                                                                                                                                                                                                          |                    |
|                           | For optimal image quality with slim patients, the shoulder coil cups should be positioned at the inside guide rails of the base plate. The outer guide rails should be used only for muscular and adipose patients. |                    |
|                           | Weight:                                                                                                                                                                                                             | 5.4 kg             |
|                           | Dimensions (L × W × H):                                                                                                                                                                                             | 489 x 444 x 46 mm  |

#### ▪ Sequences

Basic set of sequences provided on the scanner.

#### Spin echo family of sequences

SE (Spin Echo): single and multi echo

TSE (Turbo Spin Echo)

HASTE (Half-Fourier Acquisition with Single Shot Turbo Spin Echo)

SPACE (Sampling Perfection with Application-optimized Contrast using different flip angle Evolutions): for fast 3D spin echo imaging

#### Gradient echo family of sequences

FLASH (Fast Low Angle SHot): spoiled GRE (gradient echo)

VIBE (Volume Interpolated Breathhold Examination): for fast 3D gradient echo imaging

MEDIC (Multi Echo Data Image Combination)

DESS (Dual Echo Steady State)

TurboFLASH

FISP (Fast Imaging with Steady State Precession)

EPI (Echo Planar Imaging): with signal type SE or FID (Free Induction Decay)

ToF (Time-of-Flight) Angiography

PC (Phase Contrast) Angiography

CV/BEAT: comprehensive sequence for cardiac imaging

#### ▪ Computer system

| syngo MR Acquisition Workplace |                               |
|--------------------------------|-------------------------------|
| Processor                      | 1 x Quad-Core Xeon E5-1620 v3 |
| Clock rate                     | 3.5 GHz                       |
| Main memory (RAM)              | 32 GB (4 x 8 GB DDR4-2133)    |

| syngo MR Acquisition Workplace |                                                                   |                  |
|--------------------------------|-------------------------------------------------------------------|------------------|
| Hard disk drives               | 1 <sup>st</sup> hard disk (system SW)                             | 300 GB SAS       |
|                                | 2 <sup>nd</sup> hard disk (data base)                             | 300 GB SAS       |
|                                | 3 <sup>rd</sup> hard disk (images)                                | 300 GB SAS       |
| Drive controllers              | LSI 9217-4i4e                                                     |                  |
| Graphics adapter               | NVIDIA K2200                                                      |                  |
| Ethernet                       | 2 x 1 GB (1 x onboard, 1 x pcb)                                   |                  |
| Audio                          | 1 x onboard                                                       |                  |
| Optical storage                | SATA DVD burner, SATA slim DVD-ROM                                |                  |
| CD-R writer                    | Approx. 4.000 images 256 <sup>2</sup> ; DICOM Standard, ISO 9660  |                  |
| DVD-R writer                   | Approx. 25.000 images 256 <sup>2</sup> ; DICOM Standard, ISO 9660 |                  |
| Media drives                   | CD/DVD drive                                                      |                  |
| Color LCD monitor              | Screen size (diagonal)                                            | 19"              |
|                                | Horizontal frequency                                              | 31–82 kHz        |
|                                | Vertical frequency                                                | 56–75 kHz        |
|                                | Screen matrix                                                     | 1280x1024 pixels |

| MR image reconstruction system (32RX) |                                                                                                                            |
|---------------------------------------|----------------------------------------------------------------------------------------------------------------------------|
| Processor                             | 2 x Xeon E5-2609v3 (6C, 15 M Cache)                                                                                        |
| Clock rate                            | 1.9 GHz                                                                                                                    |
| Motherboard                           | X10DRH-CLN4                                                                                                                |
| Main memory (RAM)                     | 48 GB (6x8 GB DDR4-2133 ECC)                                                                                               |
| Hard disk (measurement data)          | 1xSSD 240 GB SATA, 5x300 GB SAS 10k8                                                                                       |
| Reconstruction speed                  | ≥21.000 recons per second (256 <sup>2</sup> FFT, full FOV)<br>≥96.000 recons per second (256 <sup>2</sup> FFT, 25% recFOV) |
| Parallel Scan and Recon               | Up to 12 data sets parallel + 1 high prio                                                                                  |

| PET reconstruction system (PRS)                                                                                                |                                                                                                                                                                                                            |
|--------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Provides asynchronous processing and reconstruction of acquired PET data, both at the time of acquisition and retrospectively. |                                                                                                                                                                                                            |
| Standard software                                                                                                              | Asynchronous processing/reconstruction<br>3D image-based scatter correction<br>3D iterative reconstruction<br>Fourier rebinning (FORE) + Filter-Back-Projection (FBP)<br>HI-REZ processing<br>PET Listmode |

| PET reconstruction system (PRS) |                                                                                                                                                                                           |                      |
|---------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------|
| Optional software               | Reconstruction of dynamic data sets<br>Reconstruction of cardiac gated data sets<br>Reconstruction of respiratory gated data sets<br>Cardiac PET<br>Dynamic PET<br>Respiratory PET Gating |                      |
| Processors                      | Intel Xeon E5-2630 v2 2.6 1600 6C 1 <sup>st</sup> CPU<br>Intel Xeon E5-2630 v2 2.6 1600 6C 2 <sup>nd</sup> CPU                                                                            |                      |
| Clock rate                      | 2.6 GHz                                                                                                                                                                                   |                      |
| Main memory (RAM)               | 64 GB DDR3-1866 (8x8 GB)                                                                                                                                                                  |                      |
| Video card                      | NVIDIA Quadro K600 1GB GFX                                                                                                                                                                |                      |
| Graphical processing unit (24)  | NVIDIA Tesla K40                                                                                                                                                                          |                      |
| Hard disk drives                | 1st hard disk                                                                                                                                                                             | 500 GB 7200 RPM SATA |
|                                 | 2nd hard disk                                                                                                                                                                             | 1 TB 7200 RPM SATA   |
| Media drives                    | 16 x DVD/RW drives                                                                                                                                                                        |                      |

| Acquisition control system (ACS) |                                                                                                                |                      |
|----------------------------------|----------------------------------------------------------------------------------------------------------------|----------------------|
| Processors                       | Intel Xeon E5-2630 v2 2.6 1600 6C 1 <sup>st</sup> CPU<br>Intel Xeon E5-2630 v2 2.6 1600 6C 2 <sup>nd</sup> CPU |                      |
| Clock rate                       | 2.6 GHz                                                                                                        |                      |
| Main memory (RAM)                | 64 GB DDR3-1600 (8x8 GB)                                                                                       |                      |
| Video card                       | NVIDIA Quadro K600 1 GB GFX                                                                                    |                      |
| Hard disk drives                 | 1 <sup>st</sup> hard disk                                                                                      | 500 GB 7200 RPM SATA |
|                                  | 2 <sup>nd</sup> hard disk                                                                                      | 1 TB 7200 RPM SATA   |
| Media drives                     | 16 x DVD/RW drives                                                                                             |                      |

| syngo MR Workplace (optional)                                          |  |
|------------------------------------------------------------------------|--|
| Color LCD monitor and host computer as for syngo Acquisition Workplace |  |

| Euro Connector (System IEC320) |             |
|--------------------------------|-------------|
| Load rating                    | 100–240 Vac |
|                                | 50–60 Hz    |
|                                | max. 10 A   |

## PMU Wireless Physio Control; PMU Sensor Kit

The Physiological Measurement Unit (PMU) lets you control MR measurement sequences using patient's physiological signals (ECG, respiration and pulse).

The PMU consists of the following components:

- PERU (Physiologic ECG and Respiratory Unit): ECG and respiratory sensor
- PPU (Peripheral Pulse Unit): Pulse sensor
- External trigger input/output

The following labels can be found on the PPU and the PERU components:

| PMU Technical Data         |                                           |
|----------------------------|-------------------------------------------|
| Frequency range            | 2402 MHz to 2480 MHz                      |
| Modulation parameters      | GFSK                                      |
| IEC 60601-1 classification | Class III Internally powered ME equipment |
| Applied part               | Type BF                                   |
| Ratings, Voltage:          | internally powered DC 3.7 V               |
| IP                         | No classification used: IPX0              |
| Mode of operation          | Continuous                                |
| MR compatibility           | MR-safe                                   |

| PMU Technical Data                 |                                                                                                |
|------------------------------------|------------------------------------------------------------------------------------------------|
| <b>Respiratory cushion</b>         |                                                                                                |
| Weight                             | 0.32 kg                                                                                        |
| Dimensions                         | 205 x 55 x 13 mm (excl. pneumatic tube to connect sensor/PERU)                                 |
| <b>Application cushion</b>         |                                                                                                |
| Weight                             | 0.4 kg                                                                                         |
| Dimensions                         | 400 mm x 105 mm x 21 mm                                                                        |
| <b>Respiratory belt</b>            |                                                                                                |
| Weight                             | 0.065 kg                                                                                       |
| Dimensions                         | 60x60x105mm                                                                                    |
| <b>Peripheral Pulse Unit (PPU)</b> |                                                                                                |
| Weight                             | 0.6 kg                                                                                         |
| Dimensions                         | Base unit (electronics): 138 mm x 57 mm x 19 mm<br>Attached flexible cable/tube: length 475 mm |
| <b>Charging station</b>            |                                                                                                |
| Weight                             | 0.604 kg                                                                                       |
| Dimensions                         | 138 mm x 123 mm x 35 mm (excl. USB power cable)                                                |
| Input supply                       | 5V                                                                                             |
| Charging current                   | Max. 1000 mA                                                                                   |
| Power                              | Max. 5 W                                                                                       |
| <b>LITHIUM-BATTERIE</b>            |                                                                                                |
| Weight                             | 0,026 Kg                                                                                       |
| Dimensions (mm)                    | 61 × 37.5 × 5.95 mm                                                                            |
| Nominal Voltage                    | 3.7 V                                                                                          |
| Nominal Capacity                   | 1320 mAh                                                                                       |
| Minimum Capacity                   | 1260 mAh                                                                                       |

▪ **Network/data coupling**

The MR system provides the level of safety according IEC 60950-1 outside patient environment. All equipment connected to the system's network/data couplings must also provide minimum level of safety according IEC 60950-1.

▪ **syngo MR image viewing and filming for MR and PET**

**Image recall**

Images are stored in DICOM format, allowing fast image access and recall

Combines images in a new series, e.g., post-processed images

Loads images automatically

Inline Display

**Image display**

Various display layouts selectable

Up to 3 patients can be simultaneously active in the viewer

Image annotation and labeling

Non-interpolated display

Fast paging through up to 500 images with 15 images/s for full screen display

**Windowing**

Freely selectable window width and center

Windowing on succeeding images

Auto-windowing for optimized contrast

Saves and sends window values

**Interactive movie/Automatic Movie for cine display**

Paging by dragging the mouse or Automatic Movie mode by clicking the icon

**Evaluation**

Parallel evaluation of up to 40 regions of interest

Circle

Rectangle

Freehand ROI

Pixel lens with position marker

Statistical evaluation

Area

Standard deviation

Mean value

Min/max values

Image scrolling

Magnification

Distance

Angle

**2D Post-processing**

Image manipulations

Reversal of gray-scale values

Image rotation by 90° or by user-defined angle

Flip horizontally/vertically

Image zoom and pan

Shutter

Annotation

Reconstruction of user-defined ranges of parallel, radial or freehand cuts  
Selectable slice thickness and slice increment of reconstructed images  
Storing of post-processing protocols  
Annotations and 2D evaluations such as distance and ROI

## MIP – Maximum Intensity Projection for MR images

3D reconstructions of vessels from a 3D data set, or a 2D sequential slice data set (acquired with dedicated MR Angiography sequences)

Volume of Interest (VOI) defined to increase reconstruction speed and to improve image quality  
Freehand MIP: virtually eliminating interfering overlays by manually drawn, irregular contours in all dimensions

Arbitrary views along any direction can be defined interactively with mouse-driven virtual trackball  
Multiple view angles around any orthogonal axis  
Projections displayed as single images, as interactive movie or by fast paging  
MIP thin/MIP thick

## MinIP – Minimum Intensity Projection for MR images

Similar to MIP but reconstructs the minimum intensity (e.g., for Dark Blood techniques)

## SSD – Shaded Surface Display for MR images

Three-dimensional display of surfaces, such as vessels

Selectable variable threshold values  
Multiple view angles around any orthogonal axis  
Rectangular and irregular Volumes of Interest (VOI) can be defined to improve image quality

## Fusion

Image fusion of multiple 3D data sets with alpha blending, i.e. overlay of two images with manual setting of opacity

Multiple 3D data sets from different modalities (MR, CT, Nuclear Medicine, PET)  
Visual alignment, automatic registration, or landmark based registration

## syngo Fused Vision 3D

Exceptional visualization tool for fused images. Enables 3D visualization of fused anatomical and functional volume data using advanced rendering techniques.

Special 3x3 layout for comprehensive view of functional (e.g., PET, SPECT), morphological (e.g., CT, MR) and fused images in one display  
Correlated rotating Maximum Intensity Projection (MIP) gives overview and single-click access to suspect regions visible in PET and SPECT volumes

### Specific absorption rates

The specific absorption rate for the RF power and the rate of change of the gradient fields are checked according to the requirements of norm IEC 60601-2-33.

### Cooling system

|                           |                          |
|---------------------------|--------------------------|
| Water consumption         | >100 l/min <sup>1)</sup> |
| Heat dissipation to water | ≥ 70 kW                  |

1) Water temperature: 6–12 °C (43–54 °F)

#### Line power supply

| Variable         | Values                                   | Tolerance |
|------------------|------------------------------------------|-----------|
| Voltage          | 380 V, 400 V, 420 V, 440 V, 460 V, 480 V | ±10%      |
| Frequency        | 50/60 Hz                                 | ±1 Hz     |
| Connection value | 95 kVA                                   |           |

For all products, line power has to be supplied via an on-site system contact or via another multipole shut-down mechanism. Room installation has to be in compliance with VDE 0100-710<sup>1)</sup>.

1) In all countries, compliance with local and national legal regulations is required.

However, we strongly recommend compliance with the regulations described herein – to the extent permitted by relevant local and national laws – in order to ensure the safety of operating personnel, patients, and third parties.

#### Power rating plates

**380/400/420/440/460/480 V ~**  
**3PHASE 50/60 Hz**  
**SHORT – TIME kVA: 105**  
**LONG – TIME kVA: 35**

#### Power consumption

|                                                                                                                                              |          |
|----------------------------------------------------------------------------------------------------------------------------------------------|----------|
| System off                                                                                                                                   | 8.33 kW  |
| Stand-by                                                                                                                                     | 10.18 kW |
| Ready for measurement <sup>1</sup>                                                                                                           | 22.68 kW |
| Typical examination <sup>1</sup>                                                                                                             | 32.37 kW |
| Remarks:<br>All values are typical values, applicable for 400V / 50Hz.<br>Consumption for compressor, cold head and separator pump included. |          |

1) The power consumption described herein is based on results that were achieved in a setting according to the COCIR methodology MRI - Measurement of the energy consumption (<http://www.cocir.org/site/index.php?id=46>). Since many variables impact power consumption (e.g. sequences used for scanning and sequence parameters, scan time), there can be no guarantee that each customer will achieve the same values.

#### Dimensions

| Component                                                                                                | Width [cm] | Depth [cm] | Height [cm] | Weight [kg] | Heat dissipation [kW] |
|----------------------------------------------------------------------------------------------------------|------------|------------|-------------|-------------|-----------------------|
| Examination Room                                                                                         |            |            |             |             |                       |
| Magnet 3 Tesla AS (incl. Helium)                                                                         | 215        | 200        | 222         | 6300        | n/a                   |
| Biograph mMR in operation, incl. gradient coil, body coil, PET detector, DEAs, patient table, and covers | 230        | 335        | 242         | 9000        | 2                     |

| Component                                                                                       | Width [cm] | Depth [cm] | Height [cm] | Weight [kg] | Heat dissipation [kW] |
|-------------------------------------------------------------------------------------------------|------------|------------|-------------|-------------|-----------------------|
| Patient table                                                                                   | 54 (top)   | 268        | 49-92       | 250         | 0.2                   |
| Required min. Room height clearance                                                             |            |            | 273         |             |                       |
| Min. transport dimensions                                                                       | 2301       | 3651       | 245         |             |                       |
| Control Room                                                                                    |            |            |             |             |                       |
| syngo Acquisition Workplace (table+monitor)                                                     | 120        | 80         | 117 (72+45) |             |                       |
| Host computer                                                                                   | 22         | 46         | 47          |             |                       |
| Equipment Room                                                                                  |            |            |             |             |                       |
| Electronics cabinet MR, incl. system control, RF system, gradient power system, image processor | 160        | 65         | 1982        | 1500        | ≤10                   |
| Cooling system                                                                                  | 65         | 65         | 1892        | 500         |                       |
| Electronics cabinet PET, incl. system control, image processor                                  | 160        | 65         | 1982        | 800         | <2                    |

1) Without covers and electronics

2) Without attachments

■ **Laser light localizer**

|                                                        |     |
|--------------------------------------------------------|-----|
| Laser class:                                           | 2M  |
| Maximum output of laser radiation (W/m <sup>2</sup> ): | ≤25 |
| Pulse pattern                                          | CW  |
| Wavelength (nm)                                        | 640 |

## 5.4 Labeling information

### 5.4.1 At the electronics cabinet

Name plate label

| SIEMENS                                                                                                                                                                                                                    |                                                                                                                                                                                                           |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Biograph mMR</b><br><br>(01) 00000123456789 GTIN<br>(240) 01234567 MODEL<br>(422) 012 ORIGIN DE<br>(21) 012345 SERIAL<br>(20) 01 REV |  <br>2015 0123<br><b>SYSTEM IVK</b> |
|  Siemens Healthcare GmbH, Henkestr. 127<br>91052 Erlangen, Germany<br><b>Made in Germany</b>                                            |                                                                                                                                                                                                           |

## 5.4.2 Country-specific

Each unit of the Biograph mMR delivered to Russia, Belarus, Uzbekistan and Kazakhstan will be endorsed with the EAC label.

### Artwork of the MD labeling in Russian

Additionally a Russian label will be applied with the following content after license is granted by the Russian authorities:

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### Labels according to EN 50419:2006

Within the EU, products identified with this symbol are subject to guidelines 2002/96/EG for old electrical or electronics system and 2012/19/EU respectively. Please contact Siemens Service in case of questions about returning and disposing the MR system and/or its components and accessories.


### Revision label

| Revision: |    |    |    |    |    |    |    |    |    |    |    |    |    | Model - No.: |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
|-----------|----|----|----|----|----|----|----|----|----|----|----|----|----|--------------|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| 00        | 01 | 02 | 03 | 04 | 05 | 06 | 07 | 08 | 09 | 10 | 11 | 12 | 13 | 14           | 15 | 16 | 17 | 18 | 19 | 20 | 21 | 22 | 23 | 24 | 25 | 26 | 27 | 28 | 29 | 30 | 31 | 32 |
| 33        | 34 | 35 | 36 | 37 | 38 | 39 | 40 | 41 | 42 | 43 | 44 | 45 | 46 | 47           | 48 | 49 | 50 | 51 | 52 | 53 | 54 | 55 | 56 | 57 | 58 | 59 | 60 | 61 | 62 | 63 | 64 | 65 |
| 66        | 67 | 68 | 69 | 70 | 71 | 72 | 73 | 74 | 75 | 76 | 77 | 78 | 79 | 80           | 81 | 82 | 83 | 84 | 85 | 86 | 87 | 88 | 89 | 90 | 91 | 92 | 93 | 94 | 95 | 96 | 97 | 98 |

### 5.4.3 On the magnet cover

Name plate label

#### 5.4.4 At the magnet bore (front)

Observe operator manual


#### 5.4.5 EAC Label

Each packaging unit delivered to Russia, Belarus, Uzbekistan and Kazakhstan will be endorsed with the EAC label right beside the system labels.

#### 5.4.6 Description of the symbols on the labels

|                                                                                     |                    |
|-------------------------------------------------------------------------------------|--------------------|
| GTIN                                                                                | GTIN number        |
| Model                                                                               | Model              |
|  | Serial number      |
|  | Manufacturer       |
|  | REF number         |
|  | Batch code         |
|  | Manufacturing date |

|                                                                                     |                                                                         |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------|
|    | Do not dispose this item. For disposing, please follow laws/regulations |
|    | EC conformity                                                           |
|    | EC representative                                                       |
|    | Attention! Please refer to Instruction for use                          |
|    | Please refer to instruction for use                                     |
|   | Observe operator manual                                                 |
|  | Temperature range                                                       |
|  | Humidity range                                                          |
|  | Applied part type B                                                     |
|  | Applied part type BF                                                    |
|  | Class II                                                                |

#### 5.4.7 For laser light localizer

Safety certificate label (International)

## 6 Sterility

MR system device is not intended to be sterile.

## 7 The list of International Standards and Directives applied to the Medical Device

| Standard                                                              | Description                                                                                                                                                                           |
|-----------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| IEC 60601-1:2012<br>EN 60601-1:2006 + Cor.:2012<br>+ A1:2013          | Medical electrical equipment - Part 1: General requirements for basic safety and essential performance                                                                                |
| IEC 60601-2-33:2010/A1:2013<br>/A2:2015<br>EN 60601-2-33:2010/A2:2015 | Medical electrical equipment - Part 2-33: Particular requirements for the basic safety and essential performance of magnetic resonance equipment for medical diagnosis                |
| IEC 60601-1-2:2014<br>EN 60601-1-2:2015                               | Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances - Requirements and tests |
| IEC 60601-1-6:2010<br>EN 60601-1-6:2010                               | Medical electrical equipment –<br>Part 1-6: General requirements for safety –<br>Collateral standard: Usability                                                                       |

| Standard                                                  | Description                                                                                                                              |
|-----------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| IEC 60825-1:2014<br>EN 60825-1:2014                       | Safety of laser products - Part 1: Equipment classification, requirements and user's guide                                               |
| IEC 62304:2006+A1:2015<br>EN 62304:2006+Cor.:2008+A1:2015 | Medical device software - Software lifecycle processes                                                                                   |
| IEC 62366:2007<br>EN 62366:2008                           | Medical devices - Application of usability engineering to medical devices                                                                |
| EN 1041:2008/A1:2013                                      | Information supplied by the manufacturer of medical devices                                                                              |
| EN ISO 15223-1:2016                                       | Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements |
| ISO 10993-1:2009<br>EN ISO 10993-1:2009+AC2010            | Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process                               |
| ISO 14971:2007<br>EN ISO 14971:2012                       | Medical devices - Application of risk management to medical devices                                                                      |
| ISO 13485:2016<br>EN ISO 13485:2016                       | Medical devices - Quality management systems - Requirements for regulatory purposes                                                      |

Medical

asic safe

for the




r basic  
magnet

## 8 The list of equipment and measurement tools necessary for quality control of the medical device

The following accessories for quality assurance are available at the customer. The phantoms are not calibrated.

### 8.1 Measurement phantoms

|                                                                                     |                                                                                                                                                                                                                                      |
|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | <p>Plastic bottle, 1900 ml</p> <ul style="list-style-type: none"> <li>■ Solutes per 1000 g H<sub>2</sub>O dist.:<br/>3.75 g NiSO<sub>4</sub> * 6 H<sub>2</sub>O + 5 g NaCl</li> </ul>                                                |
|   | <p>Spherical phantom D170</p> <ul style="list-style-type: none"> <li>■ Diameter: 170 mm</li> <li>■ Contents: 2570 ml</li> <li>■ Solutes per 1000 g H<sub>2</sub>O dist.:<br/>1.25 g NiSO<sub>4</sub> * 6 H<sub>2</sub>O</li> </ul>   |
|  | <p>Blue spherical phantom</p> <ul style="list-style-type: none"> <li>■ Diameter: 240 mm</li> <li>■ Contents: 7300 ml</li> <li>■ Solutes per 1000 g H<sub>2</sub>O dist.:<br/>BAYOL-Oil + 0.011 g MACROLEX blue</li> </ul>            |
|  | <p>Blue body loader</p> <ul style="list-style-type: none"> <li>■ Contents: 11000 ml</li> <li>■ Solutes per 1000 g H<sub>2</sub>O dist.:<br/>2.6 g MnCl<sub>2</sub> * 4 H<sub>2</sub>O + 2.0 g NaCl + 0.011 g SICOVIT blue</li> </ul> |

8.2

Phar

Phan

Sphe


9 7

T  
h  
T  
re  
T  
ir

9.1

•

## 8.2 Phantom holders

## 9 Transportation requirements and conditions

The magnet can be shipped and transported worldwide with or without liquid helium. Maximum possible helium level during shipping: 95%.

The magnet is classified as "Dangerous goods" as defined by the ADR, IMDG and IATA codes and their requirements. Liquid helium is classified under UN1963 (Helium, refrigerated liquid).

Transportation of the system containing liquid helium must be in accordance with the appropriate national and international codes.

### 9.1 Mechanical transport conditions

- Using an air-ride truck with or without a shipping pallet the following maximum peak acceleration measured at the magnet mounting points during transportation shall not exceed:
  - 5.0g vertical shock load applied to the cryostat mounting points.
  - 2.0g horizontal shock load applied to the cryostat mounting points in any direction.
  - Shock is defined as a dynamic force with duration of no more than 0.5 seconds.

## 9.2 Road Transport

- Road transport is allowed up to 15.000 miles.
- Road transport may be with or without an approved shipping pallet. With NO shipping pallet, the magnet must be supported using transport beams.
- Road vehicles must be fitted with air ride suspension, whether or not a shipping pallet is used.
- Always transport a magnet using vehicles fitted with air ride suspension.
- Road vehicles must be suitable for the intended usage, be properly maintained at all times, and comply with all local laws and regulations.
- Relevant requirements of the ADR code shall be satisfied.
- The vehicle and operator must be approved by Siemens.
- This specification does not constitute acceptance for use as part of a mobile system. Mobile magnets require a specific trailer design for continuous road transportation.

## 9.3 Air transport

- Air transport is allowed up to 30.000 miles.
- If air transport is used for any section of the route, an approved shipping pallet must be used for the entire transportation process.
- Relevant requirements of the IATA code must be satisfied.

## 9.4 Sea shipment

- Sea transport is allowed up to 30.000 miles.
- If sea shipment is used for any section of the route, an approved (magnet system related) shipping pallet must be used for the entire transportation process.
- Relevant requirements of the IMDG code must be satisfied.

## 9.5 Rail Transportation

- Rail Transportation is not allowed!
- Transportation through the channel tunnel (UK to France) is permissible if the system is mounted on an approved air ride truck. The air ride suspension may be inactive.

## 9.6 M

Acc

Peak

Ambi

Rela

Tem

Baro

Max.

1. Pay  
time

## 9.7 S

Acc

Pea

Ambi

Rela

Tem

Baro

Max

1. Pa  
tir

## 9.8 Lif

- Lif  
me

## 9.6 Magnet specifications

| Magnet Transport specifications including PET gantry |                                                                     |
|------------------------------------------------------|---------------------------------------------------------------------|
| Acceleration                                         | Amplitude: 2 m/s <sup>2</sup> (0.2g)<br>Frequency range: 2 - 100 Hz |
| Peak acceleration                                    | 5 g vertical, 2 g horizontal<br>(20 ms time period)                 |
| Ambient temperature                                  | +10 °C to +30 °C                                                    |
| Relative humidity                                    | 30% to 70% non-condensing                                           |
| Temperature gradient                                 | 4 °C/h                                                              |
| Barometric pressure                                  | 70 kPa to 106 kPa                                                   |
| Max. storage time                                    | 2 month <sup>1</sup>                                                |

1. Pay attention to the Helium level. Refilling might be necessary depending on helium level and transport time!

## 9.7 System components Specifications

| Storage specifications |                                                                     |
|------------------------|---------------------------------------------------------------------|
| Acceleration           | Amplitude: 2 m/s <sup>2</sup> (0.2g)<br>Frequency range: 2 - 100 Hz |
| Peak acceleration      | 5 g vertical, 2 g horizontal<br>(20 ms time period)                 |
| Ambient temperature    | +10 °C to +30 °C                                                    |
| Relative humidity      | 30% to 70% non-condensing                                           |
| Temperature gradient   | 10 °C/hr                                                            |
| Barometric pressure    | 70 kPa to 106 kPa                                                   |
| Max. storage time      | 2 month <sup>1</sup>                                                |

1. Pay attention to the Helium level. Refilling might be necessary depending on helium level and storage time!

## 9.8 Lifting and Handling

- Lifting the magnet must always be conducted by personnel who are trained in handling components and medical equipment of this size and weight.

- Lifting during transportation including the final movement onto the installed site must be performed using a suitable sized crane or forklift capable of lifting the full magnet weight.
- Suitable straps or chains capable of lifting the full magnet weight must be used. The straps or chains must be secured to the magnet lifting points using the four shackles supplied and fitted to the magnet.

## 10 Shelf life and storage

Siemens considers the subject devices exempt from having a shelf life due to the fact that the devices are manufactured upon order and is therefore not susceptible to degradation that would lead to functional failure. Once an order is placed and a machine is manufactured, it is transported, installed and QC and tune up tests are completed on site to ensure proper function of the MR-PET system. The MR-PET system does not spend time on the "shelf" that would pose a risk to the performance of the device or its components.

### 10.1 Storage specifications

Storage conditions based on system part localization:

| Storage specifications |                                                                     |
|------------------------|---------------------------------------------------------------------|
| Acceleration           | Amplitude: 2 m/s <sup>2</sup> (0.2g)<br>Frequency range: 2 - 100 Hz |
| Peak acceleration      | 5 g vertical, 2 g horizontal<br>(20 ms time period)                 |
| Ambient temperature    | +10°C to +30°C                                                      |
| Relative humidity      | 30% to 70% non-condensing                                           |
| Temperature gradient   | 10°C/hr                                                             |
| Barometric pressure    | 70 kPa to 106 kPa                                                   |
| Max. storage time      | 2 month <sup>1</sup>                                                |

1. Pay attention to the Helium level. Refilling might be necessary depending on helium level and storage time!

## 11 Expected service life

The Biograph mMR is designed for an expected service life of 10 years with operation under average conditions for which the system is expected to be safe for use (according to e.g. IEC 60601-1-2005).

Please be aware that the safety checks and their time intervals as specified in the System Owner Manual are mandatory and where appropriate, preventive measures have to be adopted or repairs performed, to ensure the safety of the system for this time.

If the products are operated for longer than this time, additional checks and possible repairs beyond the usual maintenance procedures may be necessary according to Siemens requirements to ensure functional integrity and operating safety.

## **12 Procedure of recycling and disposal (Medical Device recycling or disposal data)**

On disposing of the system or parts thereof, currently valid environmental legislation must be observed. Examples of environmentally relevant components are:

- Packaging materials
- Accumulators and batteries
- Transformers
- Capacitors
- Monitor picture tubes
- Phantoms

For details contact your local customer service representative or your Siemens regional office.

Waste class in accordance with SanPiN: 2.1.7.2790-10 - A.

System components hazardous to persons or the environment must be disposed of with care and in compliance with legally binding ordinances.

## **13 Usage conditions**

### **13.1 Installation requirements;**

To ensure correct installation, complete all preliminary steps specified.

For detailed information on preparing the room, transporting the system, and other installation requirements, consult the project manager.

The on-site project manager is responsible for all room and system-related preparations.

The following prerequisites must be met:

- The room must be fully prepared and cleaned.
- Power and water connections must be available.
- The internal resistance of the on-site power supply unit must correspond to the requirements.
- The transport route for the system must be defined and free of obstructions (regarding dimensions, access, load).
- Ensure that the magnet room is properly equipped.
- Ensure that the openings of the RF shield are sufficiently large.
- Ensure that there is sufficient space to install the magnet.
- Take account of narrow passages due to suspended ceilings, carriers, and columns.
- Determine the final location of the magnet in the magnet room.
- Ensure that the primary cold-water supply is ready for operation.

- Block off the unloading zone to ensure that passersby are not harmed.
- Do not block escape routes or essential access routes to the hospital.
- The rooms must be properly equipped. The RF-shield, including the opening for inserting the magnet, must be installed. The magnet room must also be equipped with sufficient lighting. The floor must be clear so that the magnet and system components can be transported.
- The four corners of the magnet must be drawn on the floor of the magnet room to determine the magnet's location. Refer to the project plan.
- The floor from the unloading zone to the final location must be covered with metal plates (at least 5 mm thick).
- The tools required for installation must be available on site (steel rollers, hydraulic lifts, etc.).

### 13.2 Additional maintenance information;

Regular inspections should be a part of the annual maintenance performed by Siemens Service under a service contract. We therefore recommend that you take out a service contract. It is also possible to take out a service contract based on system usage. Please contact your authorized local representative in this matter.

The name and address of the local representative in Russia:

Siemens Healthcare LLC

Dubininskaja Str. 96

115093 Moscow, Russia

Tel.: +7 495 737 12 52

[www.siemens-healthineers.com/ru/](http://www.siemens-healthineers.com/ru/)

*Customer's support center Siemens Healthcare LLC (24 hours):*

Phone: 8 800 550 50 22

Email: [service.ru@siemens-healthineers.com](mailto:service.ru@siemens-healthineers.com)

### 13.3 Information about calibration of the MD / requirements to calibration (if applicable).

There is no calibration required for the PET-MR systems. Nevertheless, there are different checks and tests which have to be done as described in the **System Owner Manual**.

## 14 Warranty

The system owner is responsible for the software licenses including software warranty conditions.

The specific warranty conditions are usually defined in the contractual agreement between Siemens Healthcare GmbH and the system owner.

## 15 Pro

### 15.1 Au

The MR-F  
necessary  
technical  
This user  
user must  
and rescu  
Please co  
frequency

#### 15.1.1

Term u

User/O  
person

System

MR wor

Siemen  
person

The syste  
owner/op  
employee  
Operating  
rescue pe  
(e.g. floor

## 15 Protective precautions for use

### 15.1 Authorized operating personnel

The MR-PET system must be operated according to the intended use and only by qualified persons with the necessary knowledge in accordance with country-specific regulations, e.g. physicians, trained radiological technicians or technologists, subsequent to the necessary user training.

This user training must include basics in MR technology, PET technology as well as safe handling of MR systems. The user must be familiar with potential hazard and safety guidelines the same way the user is familiar with emergency and rescue scenarios. In addition, the user has to have read and understood the contents of the operator manual. Please contact Siemens Service for more information on available training options and suggested duration and frequency of such training.

#### 15.1.1 Definitions of different persons

| Term used                         | Explanation                                                                                                                                                                                     |
|-----------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| User/Operator/Operating personnel | Person who operates the system or software, takes care of the patient or reads images.<br>Typically physicians, trained radiological technicians, or technologists.                             |
| System owner                      | Person who is responsible for the MR environment. This includes legal requirements, emergency plans, employee information and qualifications, as well as maintenance/repair.                    |
| MR worker                         | Person who works within the controlled access area or MR environment.<br>User/Operator as well as further personnel (for example, cleaning staff, facility manager, service personnel).         |
| Siemens Service/service personnel | Group of specially trained persons who are authorized by Siemens to perform certain maintenance activities.<br>References to "Siemens Service" include service personnel authorized by Siemens. |

The system owner manual and operator manual contain the most important safety aspects for which the owner/operator of the PET-MR system are responsible. These include legal requirements, emergency plans, employee information and qualifications, as well as requirements that must be met in the examination room. Operating personnel as well as personnel who are not regulars in the examination room (e.g. cleaning personnel, rescue personnel) must be informed about the special conditions in the examination room. Magnetizable devices (e.g. floor polishers, vacuum cleaners, wheelchairs, metal gurneys) must not be used in the examination room.

**16 Contact for reclamation**

In case of any questions or reclamations, the following shall be contacted:

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Herewith, I certify the foregoing signature,  
acknowledged before me, of

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Mrs. Velasco Alvarado is personally  
known.

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