ESSILOR
INTERNATIONAL
DIVISION INSTRUMENTS
81 Boulevard Jean-Baptiste Oudry
94000 CRETEIL

AFFIRM

Aurélie JARDIN

Head of Regulatory Affairs

Instruments Division

ESSILOR INTERNATIONAL

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03/11/2022

User's manual

Device for complex eye diagnostics, models WAM700 + and WAC700

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2022

INTRODUCTION

The apparatus for complex eye diagnostics, models WAM700+ and WAC 700 (hereinafter referred to as apparatus, device, device, product) is designed to determine optical distortions (aberrations) present in the human visual system.

Scope - ophthalmological offices of hospitals, clinics, health facilities, optics salons.

INDICATIONS

- Measurement of the refraction of the eye with both lower and higher order aberrations.
- Measurement of the shape of the cornea.
- Retro illumination of the eye.
- Measurement of intraocular pressure without contact with the eye to assess glaucoma.
- Photographing the eye and obtaining an image of the eye to assess corneal thickness.


- Drug or alcohol intoxication
- Mental illness accompanied by aggressive or inappropriate behavior

POSSIBLE COMPLICATIONS

Not identified

Design

WAM 700+/ WAC700 _


No.	Description
1	Optical measuring head Contains the hardware used to perform all measurements during patient diagnostic tests.
2	Switch on / off Shuts down the instrument.
3/6	Ventilation holes
4	Power cable connector Connect the power cable to this connector; Plug the other end of the cable into a standard electrical outlet.
5	Connectors for external devices See the table below for a list of available connectors.
7	casing Contains groups of computers and other electronics
8	LCD touch screen

The following connectors are available for connecting external devices:


Connector / port type	Usage
USB (4 connectors / ports)	• Connect an external hard drive or flash drive to export data from the group database and/or to import data for executables; • Connect a keyboard to use instead of or along with a virtual keyboard; • Connect a mouse to use instead of or along with the touch screen; • Connect a printer to print on paper.
Network (2 connectors / ports)	Connect the network cable to connect the device to the local network. You will be able to export and import data from computers over the local network.
Port (RS-232)	Connect a compatible optometric device such as a phoropter to the device
VGA port / connector	Connect an external monitor to the unit.

Patient side

No.	Description
1	Headrest The patient should rest their forehead on the headrest in all diagnostic procedures.
2	Chinrest The patient should rest their chin on the chin rest in all diagnostic procedures.

Print side

No.	Description
1	ventilation slot
2	LED printer Printer status indicator: • OFF indicator : The printer will be turned off. • ON indicator (not blinking): The printer is on and ready to print. • The indicator is blinking: The printer is on. There may be a problem or there is no paper.
3	Printer paper feed button Press to feed a small amount of paper through the paper slot.
4	paper slot Paper is fed through this slot when it is loaded.
5	Paper roll compartment

Pin for face stop	
	
Network cable 220V	dust cover
	
B chin rest pad	Paper for built-in thermal printer
	

SPECIFICATIONS,

Machine Specifications

	WAM700+	WAC700
Supply voltage	100-240VAC, 50/60Hz	100-240VAC, 50/60Hz
Power consumption, no more	300 W	300 W
Electric current protection type	Class 1	Class 1
Degree of protection against electric current	Type BF	Type BF
IP classification	IP X0	IP X0
Size (W × D × H)	312×530×570±2mm	312×530×570±2mm
W8	25±0.2 kg	25±0.2 kg

Operating conditions

	WAM700+	WAC700
temperature	+ 10 ° C to + 35 ° C	+ 10 ° C to + 35 ° C
hygrometry	30% to 90%	30% to 90%

Storage conditions

	WAM700+	WAC700
temperature	-10°C to +55°C	-10°C to +55°C
hygrometry	10% to 95%	10% to 95%

Transport conditions

	WAM700+	WAC700
temperature	-40°C to +70°C	-40°C to +70°C
hygrometry	10% to 95%	10% to 95%

Specifications

	WAM700+	WAC700
Spherical power range	-20D to +20D at 12mm tip distance	-20D to +20D at 12mm tip distance
Cylinder power range	0D to 8D	0D to 8D
Axis	0 - 180°	0 - 180°
Measurement area	2.0 mm - 7.0 mm pupil (3 zls)	2.0 mm - 7.0 mm pupil (3 zls)
Working distance	94 mm	94 mm
Method	Shaka -Hartmann	Shaka -Hartmann

Topography of the cornea

	WAM700+	WAC700
Number of rings	24	24
Number of measuring points	6144	6144
Number of analyzed points, not less than	100000	100000
Diameter of covered corneal area at 43D	0.75 mm to 10 mm	0.75 mm to 10 mm
Diopter	37.5D - 56D	37.5D - 56D
Repeatability	0.03mm	0.03mm
Method	Disc Placido	Disc Placido

Specifications of the pachymeter for the WAM 700+ model

Measurement Wavelength	455 nm without UV
measurement range	50 μm - 1300 μm
Accuracy, no more	5 μm
Method	Static Scheimpflug

Specifications of the tonometer for the WAM 700+ model

Calibration range	7 mm - 50 mmHg .
Accuracy	+/- 2 mmHg
Method	contactless tonometry

Cataract Screening Specifications

	WAM700+	WAC700
Measurement Wavelength	850 nm	850 nm
Method	retro illumination	retro illumination
Beam size	1.2 $\pm$ 0.2 mm	1.2 $\pm$ 0.2 mm
Laser output power	12 $\pm$ 2 μW	12 $\pm$ 2 μW

LCD Screen Specifications

	WAM700+	WAC700
Screen type	TFT LCD	TFT LCD
The size	10" WSVGA resolution (1024 $\times$ 600 pixels)	10" WSVGA resolution (1024 $\times$ 600 pixels)
Contrast	500	500

The product must maintain its performance under vibration loads with the following parameters: frequency range 10-55 Hz with a displacement amplitude of 0.15 mm.

The product, during transportation, must be resistant to mechanical stress, with the following parameters:

- vibration loads - frequency range 10-55 Hz with a displacement amplitude of 0.35 mm;
- impact loads - peak impact acceleration 100 ms⁻² with duration of impact acceleration 16 ms , impact frequency up to 120 min⁻¹.

The time for establishing the operating mode should be no more than 1 min.

The mode of operation must be long.

The product must be classified as a class I ophthalmic device.

The product must use fuses with a rating of T2A L 250V.

The external surfaces of the parts of the apparatus must be resistant to repeated chemical disinfection using a 3% hydrogen peroxide solution with the addition of a detergent.

The average service life should be at least 5 years. The criterion for the limit state of the product is the impossibility or inexpediency of restoration.

The shelf life of the product from the date of manufacture should be 3 years.

Overall dimensions and w8 of comp1nts:

- Network cable 220V - cord length - 1800 $\pm$ 5 mm
- Dust cover - (W $\times$ D $\times$ H) 350 $\times$ 570 $\times$ 640 $\pm$ 2 mm
- Paper for chin rest - 120 $\times$ 35 $\pm$ 1 mm
- Pin for face stop - 20 $\times$ 3 $\pm$ 0.1mm
- Paper for the built-in thermal printer - 60 $\times$ 55 $\pm$ 1 mm

Mass of comp1nts:

- Network cable 220V - 270 $\pm$ 10 g
- Dustproof cover - 160 $\pm$ 4 g
- Paper for chin rest - 120 $\pm$ 3 g

- Pin for the face stop - 3 ± 0.2 g.
- Paper for the built-in thermal printer - 86 ± 2 g

Software:

Model	Name of the software	Identification name of the software	Version number (identification number) of the software	Digital software identifier (checksum of the executable code)
WAM700+	VX120	VX120.exe	2.0.1510.801	N/A WinXP exe
WAC700	VX110	VX110.exe	2.0.1510.801	N/A WinXP exe

Safety software class - Class A

Precautionary measures

Important!

- Do not place or use the machine in direct sunlight.
- Do not expose the machine to excessive dust or humidity.
- Do not place the machine in hot air (eg over a heater).
- Do not block the ventilation holes.
- Never place the machine close to the following types of equipment, which may interfere with the reception of commands from the remote control:
 - o Halogen lamp (direct or indirect)
 - o Fluocompact lamp
- Keep the screen surface clean. Protect it from dust, fingerprints and shock.
- When you turn off the unit, wait at least 5 seconds before turning it on again.

Procedure for operating the instrument.

Installation Procedures

Best results will be obtained if the device is placed in a room with limited lighting.

Unboxing

To unpack the instrument:

- Remove the straps around the box.
- Open the box with care.
- Lift the protective cover, remove the accessories listed in the packing list.
- **Important!**
Do not hold the screen while taking the device out of the box.
- Remove the machine from the box and place it on a table.
- Remove the protective plastic cover to open the device.

Loading paper in the printer

If roll paper is being used, the printer's LED blinks.

To insert a paper roll into the printer:

- Lift the handle at the center of the paper compartment and pull the cover down.
- If there is an empty roll left in the printer, remove it.
- Insert a new roll.
- Press the paper feed button at the top of the compartment.

Block on/off

To turn on the device:

- Press the On / Off switch.

To turn off the device:


- Using the software, in the Home tab Screen ", press the Off button .
-OR-
- Via the software, in the " Configuration " tab screen ", click the Off button .
-OR-
- On the block, press the On / Off switch .

Regardless of which of these methods you use, the device will turn off.

Software Overview

home screen


The **home** screen is the main menu of the machine. This gives you access to all interface features.

No.	Description
1	Add Patient Button Opens the patient screen with patient information
2	Measurement button Opens the measurement screen
3	Open Patient Screen Button Opens the patient screen and displays a list of patients
4	Off button Shuts down the instrument
5	Screen adjustment button Opens the Configuration Screen
6	Screen saver button Changes the splash screen on the screen.

Patient screen


The patient screen gives you access to the patient database. You can create new patient records, modify existing information, view saved results, and initiate new diagnostic procedures for a selected patient from this screen.

No.	Description
1	home button Opens the home page.
2	Patient List List of patients in the database.
3	Last name filter The list of patients is filtered to display only patients whose last names begin with the specified letter sequence.
4	Mark measurements Checking all displayed measurements.
5	Name filter The list of patients is filtered to display only patients whose names begin with the specified letter sequence.
6	List of measurements List of test results for the selected patient.
7	Clear filter button Clears all patient filters and displays all patient records in the database.
8	Mark all patients Check the checkboxes of all displayed patient records.
9	Displayed patient record Blue indicates that a patient record is selected.
10	Flagging Displays the selected patient record.
11	Screen saver button Change the screen saver.
12	Settings button Opens the configuration screen.
13	Add a Patient Button Displays information about the patient.
14	Edit Patient Data Button Opens patient data.
15	Delete Patient button

	Deletes all patient records that are checked.
16	Cut button Removes all measurements that are flagged from the selected patient record and saves them to the clipboard so they can be pasted into another patient record.
17	Delete measurements button Deletes all dimensions that are checked.
18	Export button Exports the selected measurements to an external file and prints the measurement results.
19	Download button Opens the selected measurements in the configuration screen.
20	New diagnostics button Initiates new diagnostic procedures for the selected patient. Diagnostic results are automatically stored in the patient's record.

When exporting (by clicking on the No. 18 "Export" button), the following pop-up window will appear from the patient screen:

No.	Description
1	Data transfer Options for transferring data from the device to other devices.
2	External device Export data from the machine to another device.
3	Ticket Allows you to print the data on the ticket in the printer's internal device.
4	Screen printing Check the box to print a summary screen shot, this requires connection to an external printer.
5	Send Confirms the selection and performs the transfer.


The patient screen export pop-up allows you to export data via an external device. Check the box next to the external device, and then click SEND. To print a screenshot of the screen you are viewing, select Print Screen.

Screen "Measurement"

The Measurement screen is used to select, initiate, and track the results of diagnostic tests.

Before Measurement:


When you first open the "Measurement" screen, before running diagnostics, the screen looks like this:

No.	Description
1	Home button Opens the main screen.
2	A patient The name (or ID) of the current patient. Note: If no patient has been selected before opening the measurement screen, a temporary name is generated by the system by concatenating "Temp " with a number (eg "Temp 1184").
3	Up arrow Raise the head of the device so that it is better aligned with the patient's eyes.
4	Message field A text field in which the system displays messages to you, such as instructions for you to provide a patient or information about the testing process.
5	Left eye indicator When active, indicates that the head is aligned with the patient's left eye Note: In the picture above, the right eye indicator is active (see #8) while the left is inactive.
6	Diagnostics available List of available diagnostic tests. Drag up and down the list to scroll if necessary and select the type of test you want to run.
7	Arrow input

	XML database . Just click the blue arrow to download patient data from your database. If the arrow is greyed out, it means that the input folder is empty and there is no information available for download.
8	Eye selector Select an eye or eyes to test.
9	Right eye indicator When active, indicates that the head is aligned with the patient's right eye. Note: In the picture above, the indicator of the right eye is active, while the left is not (see #5)
10	Camera view Shows the image seen through the camera. Select anywhere in the image to move the device head to align that spot with the crosshair in the center of the image.
11	Left Arrow Move the device head to the left to better match the patient's eyes.
12	Job title Select " Adult " to move the head and chin to the default adult positions, or select " Child " to move them to the default child positions. Once you've done that, you can fine-tune your head and chin position with the chin controls (#12), arrows (#3, #10, #16, and #17).
13	Chinrest Raise or lower your chin.
14	Options Open a dialog where you can change the number of WF and Tono tests to run during this diagnostic. Note: If the WF and Tono tests are not included in this diagnostic, changing the values in the dialog box will have no effect.
15	Cancel Cancel the measurement procedure and return to the previous screen.
16	To begin Start a diagnostic test.
17	Arrow to down Lower the head of the device so that it is better aligned with the patient's eyes.
18	Right Arrow Move the device head to the right to better match the patient's eyes.
19	Measurements to be taken Specifies which measurement procedures will be performed for the diagnostic selected in the Diagnostics list (No. 6)


Clicking on the option (No. 13) displays the following dialog box:

No.	Description
1	Demo version Select to demonstrate tonometer air puff
2	WF measurements Number of wavefront measurements to perform
3	Tono measurements (for WAM 700+ model) Select to change the number of tono measurements to perform
4	Reading Distance Close view distance, 30 to 60 cm
5	Cancel Closes the dialog without saving
6	Save Saves user settings

During measurement


After you select "Go" to start the selected measurement procedure, the "Measure" screen displays information about the measurement process. During the measurement process, you can adjust the position of the head and chin, just as before the start of the process. You can also speed up the measurement process, skip parts of the test, or stop the test.

No.	Description
1	Home button [Button disabled when diagnostics in progress]
2	A patient The name (or ID) of the current patient. Note: If no patient has been selected before opening the measurement screen, a temporary name is generated by the system by concatenating " Temp " with a number (eg " Temp 1184").
3	Up arrow Raise the head of the device so that it is better aligned with the patient's eyes.
4	Message field A text field in which the system displays messages to you, such as instructions for you to provide the patient or information about the measurement process.
5	Left eye indicator When active, indicates that a measurement is currently being taken on the patient's left eye. Note: In the above picture, the right eye indicator is active (see #9) while the left is not active.
6	Sensor type During WF measurement, the Shack -Hartmann image is displayed in this area.
7	Measurement in progress Indicates that this measurement is currently in progress.
8	Completed test Indicates that this measurement has already been taken.
9	Right eye indicator When active, indicates that a measurement is currently being taken on the patient's right eye. Note: In the picture above, the right eye indicator is active, while the left is not (see #5).
10	Camera view Shows the image that is displayed through the camera.

	Note: To manually center the image, select a point that should be in the center of the image. The block head moves so that this spot is aligned with the crosshairs in the center of the image.
eleven	Left Arrow Move the device head to the left to better match the patient's eyes.
12	Position Select " Adult " to move the head and chin to the default adult positions, or select " Child " to move them to the default child positions. By default, PD is also set to adults or children when you select 1 of these options. Once you've d1that, you can fine-tune your head and chin position with the chin controls (#13) and arrows (#3, #11, #17, and #18).
13	Chin Select the arrows to raise or lower your chin.
4teen	Quick Mode Complete the process of measuring up. This option is only available during WF measurements. This is useful if the patient has trouble sitting for an ex10ded period of time. Measurement results may be slightly less accurate than in standard mode.
fifteen	Skip Cancel the current step of the measurement procedure and proceed to the next step. For example, if a measurement is being taken in the patient's left eye, selecting this button will stop the measurement in that eye and start measuring the other eye.
16	Stop Stop the entire measurement process. When this button is selected, the measurement procedure is interrupted.
17	Arrow to down Lower the head of the device so that it fits better with the patient's eyes.
8een	Right Arrow Move the device head to the right to better match the patient's eyes
19	Measurement status Indicates which measures are included in the current measurement procedure and shows the status of those measures that have already been completed.


Immediately after the measurement, the Ringer edit screen will be displayed (only if a topography measurement has been taken). On this screen, the user can edit the rings detected by the device and those rings that were not detected correctly.

No.	Description
1	Erase To erase individual rings in the topographic map image, select "Erase". A list of ring positions will open. Select the position of the ring you want to erase and then move the cursor over it and only the ring at the selected position will be deleted. Note that the outer rings are in higher positions starting at position 24 and below, while the inner ring is at 1. Note: If you remove more than 10% of the topography rings, you may reduce the accuracy of the results.
2	Add Select the ring position number (displayed in blue), then follow the ring at that position, the device "draws" a new ring. Note that rings can only be drawn if the topo image is relatively transparent.
3	Automatically Allows you to remove rings by touching the screen. This is the default mode.
4	Cancel cancel changes
5	Repeat Remeasurement
6	Proceed Checks the changes made and calculates the topographic map

Call editing

If you have selected ADD or ERASE in the menu above, the following screen will appear:

No.	Description
1	Control buttons Select these buttons to add, remove, cancel, or continue.
2	Ring positions Select the position of the ring you want to change. Once selected, the ring will be blue.
3	Chosen ring Once you have selected the Ring (#2) it will be blue and can be erased/added. To erase a ring, move the cursor over the areas you want to erase. To add to the ring, move the cursor over the area you want to draw. In the image above, the selected bell is selected in position #11.


To erase the rings at will, select the AUTO mode and then slowly move the cursor over the rings you want to delete. This action can erase several rings at once. After changing the rings, select "Continue" and the changes will be saved.

Note:

Only changes made immediately after verification will be saved.

Repeat measurement

After the results screen is displayed, you can repeat the measurement by pressing the "Measure" button (see Results Screen). After pressing the "Measurement" button, the following dialog box opens:

No.	Description
1	Choose all Select all measurements for both eyes
2	Right Measurements for the right eye
3	left Measurements for the left eye
4	Both eyes Measurements for both eyes
5	List of measurements List of measurements that have been taken before and can be repeated
6	Confirm Checks the selection and performs a new measurement according to the selection

Select the measurements you want to repeat and continue the measurement as usual.


Results Screen

The Results screen displays the results of the diagnostic measurements. It opens when a diagnostic test is completed or when a test result is saved from the patient screen.

The Results screen first displays a summary which gives a condensed view of the results. The other tabs on the Results screen show some aspects of the results in more detail.

Results screen overview

The Results screen is dynamic and can have different layouts and data according to the selected diagnostic. This is a general overview of the screen of possible results:


No.	Description
1	Home button Opens the main screen.
2	Patient name, date and time Surname and name of the patient, date and time of measurement
3	Left eye indicator The results displayed on this side of the screen are those of the left eye.
4	VD Apex distance, this button displays the refraction measurement distance, it allows you to switch between default and zero.
5	Button D / mm This button allows you to switch between displaying results in millimeters and diopters.
6	Measurement button Opens the measurement screen with the selected patient name so that you can perform additional tests that were not included in the current diagnosis.
7	Export button Exports the test results to a file on an external device and/or prints the test results. For more export options, click the Export from Results tab.
8	exit button Closes the results screen and displays the default screen (which is selected in the configuration settings).
9	Tabs There may be fewer tabs in this section, depending on the selected diagnostic. Select a tab to display a topic in more detail.

Dynamic Pivot Tabs

Dynamic Summary tab is displayed by default, opening the **Result screen**. This gives a compressed view of the results.


Note:

If no results are available for a particular test type, the result fields are left blank. This chapter will show some summary tabs and their descriptions.

Summary tab for all diagnostics:


No.	Description
1	PD Visual distance of the patient (if available).
2	Refraction Displays photopic and mesoscopic refractions with pupil sizes. The image shows a photopic (above) and a mesopic (below) vision simulation. Clarity E indicates how well the patient sees from the eyes under photopic or mesopic conditions. Double-click (or double-click) E to open the day/night visual acuity simulation.
3	Add Add displays myopia and distance (in cm) measured at.
4	Topography Displays a topographic map of the eye as well as keratometry values . Double click (or double click) to open the topography map sub -tab for the eye.
5	pachymetry Displays an image of the Scheimpflug cornea as well as the thickness of the cornea at its center. Double click (or double click) to open the Biometry tab for the eye.
6	Tonometry (for WAM 700+ model) Displays glaucoma screening results, which include: tonometry test results in mmHg , mean tonometry test results, and diaphragm angles and anterior chamber depth. In addition, it displays corrected tonometry results, which are based on the formulas selected in the settings menu. Double click (or double click) to open the Tonometry tab Tab
7	Right eye indicator Indicates the test results on this side of the screen for the right eye.

Pivot table:


No.	Description
1	Refraction Displays photopic and mesoscopic refractions with pupil sizes.
2	Keratometry Displays keratometry readings , KPI values (Keratoconus index) and high order and low order aberrations
3	E simulation Clarity E indicates how well the patient sees from the eyes under photopic or mesopic conditions. Double-click (or double-click) E to open the day/night visual acuity simulation.
4	Topography Displays a topographic map of the eye. Double click (or double click) to open the topography map sub -tab for the eye.
5	D / mm Selection button to switch between diopters and mm
6	Subtitles Select subtitles to view the item in more detail

Summary table for diagnosing cataract / cornea:


No.	Description
1	Refraction Displays photopic and mesoscopic refractions with pupil sizes.
2	Keratometry Displays keratometry readings , KPI (Keratoconus Index) values, and high order and low order aberrations.
3	E simulation Clarity E indicates how well the patient sees from the eyes under photopic or mesopic conditions. Double-click (or double-click) E to open the day/night visual acuity simulation.
4	Topography Displays a topographic map of the eye. Double click (or double click) to open the topography map sub -tab for the eye.
5	pachymetry Displays a cross section of the anterior chamber as well as pachymetry for the cornea
6	D / mm Selection button to switch between diopters and mm
7	Subtitles Select subtitles to view the item in more detail

Summary table for diagnosing glaucoma:


No.	Description
1	Tonometry (for WAM 700+ model) Displays IOP readings, average IOP and corrected IOP (according to the formula selected in the setting)
2	Pachymetry (for WAM 700+) Displays cross section of cornea and pachymeter values
3	Subtitles Select to view the item in more detail.

Summary screen for WF


No.	Description
1	Refraction Displays photopic and mesoscopic refractions with pupil sizes.
2	E simulation Clarity E indicates how well the patient sees under photopic or mesopic conditions. Double-click (or double-click) E to open the day/night visual acuity simulation.
3	PSF Displays PSF uncorrected and corrected LOA
4	aberrations Displays values for LOA and HOA under photopic and mesoscopic conditions
5	Subtitles Select subtitles to view the item in more detail

Topo summary screen

No.	Description
1	Keratometry Displays K1 and K2 values
2	Aberration diagram Displays LOA and HOA values as well as KPI values (keratoconus)
3	Toro Displays a map of the topography of the cornea
4	K reading Displays keratometry values based on pupil diameter (3.5.7 mm)
5	Subtitles Select subtitles to view the item in more detail

Summary Screen for Pachymetry Diagnostics (for WAM 700+)

No.	Description
1	Pachy Image (pachymetry image) Image of the anterior chamber of the cornea with pachymetry results
2	Negative Pachy Image (negative pachymetry image) Negative image of the anterior chamber of the cornea with the results of pachymeters
3	Subtitles Select subtitles to view the item in more detail

Consolidated screen for refraction

Temp 3178 18-Mar-2015 17:30:22

PD: 70.5 mm

1 Refraction

Pup.	S	C	A
3.0	3.25	-3.25	175°
5.0	3.50	-3.25	174°

Pup.	S	C	A
3.0	-1.25	-0.25	126°
5.0	-2.00	-0.25	129°

3 VD: 12.0 mm

2 Topography

K1: 39.50D @172°	K1: 40.00D @132°
K2: 41.25D @82°	K2: 42.25D @42°
Cyl: -1.75D @172°	Cyl: -2.25D @132°
AVG: 40.38 D	AVG: 41.13 D


4 D/mm

Summary Exit

No.	Description
1	Refraction Displays photopic and mesoscopic refractions with pupil sizes.
2	Keratometry Displays keratometry readings , KPI values (Keratoconus index) and high order and low order aberrations
3	VD Displays the distance to the top of the results.
4	D/ mm Selection button to switch between diopters and mm
5	Subtitles Select subtitles to view the item in more detail


Summary tab for custom diagnostics

Custom diagnostics can contain many different measurements according to your preferences. If Topo , Pachi and Tono are selected , the resulting screen will look like this:


No.	Description
1	Topography Displays a topographic map of the eye as well as keratometry values . Double click (or double click) to open the topography map subtab for the eye.
2	Pachymetry (for WAM 700+) Displays an image of the Scheimpflug cornea as well as the thickness of the cornea at its center. Double click (or double click) to open the Biometry tab for the eye.
3	Tonometry (for WAM 700+ model) Displays glaucoma screening results, which include: tonometry test results in mmHg , average of tonometry test results. If the results are outside the average range, aperture ions are displayed here.
4	IOPC These are corrected tonometry results, which are based on the formulas selected in the settings menu. Double click (or double click) to open the Tonometry tab Tab .
5	Subtitles Select subtitles to view the item in more detail

Summary table for the diagnosis of pulpometry

No.	Description
1	pupil measurement Displays pupil measurement under different lighting conditions: scotopic , low, mesopic , photopic .
2	Pupil image Displays a pupil image for each eye.
3	Pupil Select to display the pupil subtitle, including a short video of the eye's reaction to changing lighting conditions.
4	Diagram This chart displays the pupil size for each eye under varying lighting conditions.
5	R or L and R/L Select 1 eye or both eyes at the same time.


If you switch the "R or L" button, the following screen appears:

No.	Description
1	Play / pause button Select to play, pause, rewind, or fast-forward the student's video.
2	R or L and R/L Select 1 eye or both eyes at the same time.
3	Pupil measurement Displays pupil measurement under different lighting conditions: scotopic , low, mesopic , photopic .
4	Diagram This chart displays the pupil size for each eye under varying lighting conditions.

Special Warnings

Summary tab will display warnings when results are outside the expected range:

No.	Description
1	Refraction Refraction is displayed in red: • When the difference in the area of refraction between the eyes is greater than 2D, or the difference in the cylinder between the eyes is greater than 2D. • If the difference in refraction of the same eye in mesopic and luminous light is different from 1D sphere and 0.75D cylinder • If the Sphere is greater than $\pm 20D$ or the Cylinder is greater than $\pm 8D$. • When the pupil diameter is less than 2.5 mm in at least 1 eye. • When the difference in pupil size between OD and OS is greater than 1 mm • When the cylinder is greater than 4D in either eye.
2	Add Add is displayed in red when the difference between the expected addition for age and the measured addition is other than 1D.
3	E simulation Simulation E is highlighted in red when: • Delta between night and day sphere $> 1D$ and/or cylinder greater than $> 0.75D$ • The ocular equivalent sphere of the TSF in the RMS area is greater than 0.5D in the right or left eye.
4	RMS RMS is displayed in red if the corneal or RMS equivalent sphere is greater than $> 4D$ KPI or when the corneal RMS equivalent sphere is greater than 0.5D in either eye.
5	KPI The KPI result is displayed in red when its value is higher than the normal population
6	IOPC The IOPc value is displayed in red, when the interocular pressure value is higher than 20mmHg, aperture angles are displayed.
7	Pachymetry (for WAM 700+)

Pachymetry is displayed in red when the thickness of the central cornea is lower than 400 μm or higher than 700 μm .


Maps Tab




Maps tab displays maps of the wavefront of aberrations and topographic maps of the cornea



- Eyepiece
- Cornea
- Interior
- Compare



Eyepiece

subtab displays ocular aberration wavefront maps for l or both eyes. Display options allow you to view results for photopic and mesopic conditions, view maps of both eyes at once or each eye separately, isolate lower and higher order aberrations, change the pitch and center value of the map.

No.	Description
1	Photo printing / Mesopic selector Select Day to display daylight maps or Night to display maps for nighttime lighting conditions.
2	The note Select this to view an existing note or add a note to a patient file. For more information, go to Adding a Patient.
3	Name and date Patient name, date and time of measurement
4	Story Displays map results from previous and current topographic maps according to the selected date. 
5	R or L mode Displays a single map, for the eyes selected in the eye selector.  R or L mode
6	RR or LL mode Displaying two maps as for the same eye, namely the eye selected in the eye selector; displaying two maps of the same eye allows two types of aberrations to be compared. For example, you can display HOA data on 1 map and LOA data on another.  RR or LL mode
7	RL Mode Displays two cards, 1 for each eye. Both maps are displayed with the same display options.  RL Mode



8	3D: Displays 1 card for the eye selected in the eye selector in 3D. Note:  In 3D mode, you can change the orientation by dragging your finger (or the cursor if a mouse is connected to the device) on the image. Z1s: Select this button to overlay z1rules (3, 5, and 7 mm) on the map display or remove them when they are displayed.  Map with z1rules
9	Eye indicator Points and switches between left and right eye
10	Map type selector Choose from the following card types: • Total: Display compositions of all aberrations • HOA: display only high-level aberrations • LOA: display only low level aberrations • Sphere: display only the aberrations of the sphere • Cylinder: Show only cylinder aberrations.
11	Date tab Opens a map from previous explorations. This is only available if the patient has previous exams saved in the database.
12	Z1s: Select this button to overlay z1markers (3, 5, and 7 mm) on the map display or remove them when they are displayed.  Map with z1s



13	Angles: Tilt angle guides on the map screen, or remove them when they are displayed.  Map with angle markers
14	Measurement Select to remeasure. Warning, if you choose to exit and save measurements for the first time, your previous results will be deleted.
15	Crosshairs: Overlay or remove crosshairs on the map when they are displayed  Map with crosshair guidelines
16	Autoscaling : Turn autoscaling on or off. When autoscaling is enabled, the pitch and center value of the map are set automatically, while the pitch and center value selectors are disabled. When it is deactivated, the step and center value of the selector are enabled.
17	Default Button Resets the pitch and center value to what is defined as the default value in the settings menu.
18	Central value selector Choose an average value for the color scale. Select + to increase the value or - to decrease it. Changes made to Central Value are implemented when you select Confirm .
19	Pitch selection Select the rate at which the color bar changes. Select + to increase the step or - to decrease it. Changes you make to the Step value are applied when Apply is selected .
20	Confirm button Implements the changes you make to Step and Central values .
21	Map Aberration map
22	Scale Shows the color coding scale used on the map




corneal



Topo tab displays topographic maps of the corneas of 1 or both eyes. Display options allow you to view maps from both eyes simultaneously or each eye individually, to view 4 different types of topographic maps, change map pitch and center value, and overlay data and recommendations of various types on maps.

No.	Description
1	The note Click to view an existing note or add a note to a patient file. For more information, go to adding a patient.
2	Name and date Patient name, date and time of measurement
3	story Displays map results from previous and current topographic maps according to the selected date. 
4	R or L: Displays 1 card, for the eyes selected in the eye selector.  R or L mode
5	RR or LL:

	Displays two maps as if for the same eye - namely, the eye selected in the eye selector; displaying two maps of the same eye allows you to view the same data in different ways. For example, you can display HOA data on 1 map and LOA data on another.  RR or LL mode
6	RL: Show two cards, 1 for each eye. Both cards are displayed the same way to show options.  RL mode
7	3D: Display 1 map, for the eyes selected in the eye selector, in 3D mode.  3D Mode Note: In 3D mode, you can change the orientation by dragging your finger (or the cursor if a mouse is connected to the device) on the image.
8	Value button Shows or hides values on maps  Values displayed on the map
9	Diopter/mm button Toggles the display of values between diopters and millimeters.
10	Transparent Displays a photo of the eye in the background with a projection of Placido's ring
11	Map Displays a topographic map of the cornea according to the map selected from the drop-down menu (in #5).
12	Rings Displays an image of Placido's rings projected onto the cornea. Choose


	Rings to view the call and make changes.  concentric rings
13	Kerato Shows / hides central keratometry values .
14	Pupil Overlay a circle directive on the map indicating where the pupil is or removes when it is displayed. 
15	Z1s Superposition guide z1s (3, 5 and 7 mm) on the map display or remove them when they are displayed. 
16	corners

	Move the corner pointer over the map screen or remove them when they are displayed. 
17	Cross Superimpose cross guidelines on the map display or remove them when they are displayed . Each gap has a size of 1 mm. 
18	Auto scale Activate or deactivate autoscaling . When autoscaling is enabled, the pitch and center value of the map is set automatically, and the pitch and center value of the selector are disabled. When it is deactivated, the step and the center value of the selector are enabled.
19	Map type selector Choose from the following card types: Axial: display an axial topographic map; Instant: displaying a tangential topographic map; H8: displaying a topographic map of the place; Refraction: Displays a refractive topographic map.
20	Date tab If more than 1 measurement has been taken on the current patient, this tab opens a drop-down menu with the dates of other selected measurements. See a card of a different date, select it from this menu.
21	Central value selector Choose an average value for the color scale. Select "+" to increase the value, or "-" to decrease it. Changes made to the Step value are implemented when Apply is selected .
22	Step selector Select the rate at which the color ramp changes: the range of values represented by each individual color. Select "+" to increase the step or "-" to decrease it. Changes made to the Step value are implemented when you select Confirm.
23	Confirm button Apply your changes to the Step and Central values value
24	Default button Reset Step and Central Value to default values.

25	Scale Specifies the color coding scale used on the map.
26	Map Detailed information about the selected map (topography, aberrometry or refraction).


Interior

The **internal** aberration map represents the differences between corneal aberrations and ocular aberrations for each eye. These aberrations can originate anywhere along the optical system from the posterior cornea to the retina; most of it from a crystalline lens. Maps can be displayed in both photopic and mesopic conditions.

No.	Description
1	Left eye indicator The results displayed on this side of the screen are those of the left eye.
2	Photo printing / Mesopic selector Select "Day" to display maps for daylight conditions or "Night" to display maps for night lighting conditions.
3	Map selection Select a map from this menu to display. Available cards for selection: Total, LOA, HOA, Cylinder and Coma .
4	Eye selector Select the eye you want to display.
5	Cornea Map Displays a corneal map that contains corneal aberrations.
6	Internal map Displays an internal map containing internal aberrations.
7	General Map Displays a map containing total aberrations for the selected eye.


Internal maps are useful for indicating the location of aberrations, for example :

The internal map here illustrates how most of the total aberrations are due to internal aberration. This could be a sign of a cataract or an internal eye problem.

Comparison

The chart comparison feature compares two charts for the same patient and generates a delta chart to illustrate change over time. Corneal and ocular maps are available for comparison.

No.	Description
1	Current Tab Date By default, the date selected here will be for the current measurement. To select a different date, click on the arrow and a list of previous measurements becomes available for selection.
2	Previous Date Tab Select the date of the map by which you want to compare the current map.
3	Map selection

	Select the map type from this menu to be displayed, or Ocular or Corneal maps are available for selection.
4	Map selection Select a map from this menu to display. The maps available for selection for Ocular maps are Total, LOA and HOA maps. Corneal maps are available for selection: Axial, tangential, summary, HOA and LOA maps.
5	Data at certain points To display data changes anywhere on the map, click any location on the map and the data will be displayed here.


Cornea Tab

The Cornea tab displays topographic and biometric data about the cornea

Topo Data

The Topo Data sub-tab in the Cornea tab displays information about the shape of the cornea.

The subtab displays the results of the Topo test in three tables: Keratoconus, Geometry, and details. The last table provides information on the eccentricity and radius of curvature of the cornea in various predefined positions, the radius of curvature at the meridians, and corneal keratometry.

No.	Description
1	Keratoconus Probability Index Screening results for keratoconus.
2	Geometry table Corneal geometry ($P = 1st^2$): P: corneal shape factor e: Corneal eccentricity
3	Detailed menu

	Select the data to display in the "Details" table: "Eccentricity", "Sagittal Radius", "Meridians", or "SIM-K"
4	Detailed table Displaying the data selected in the "Details" menu

Keratoconus table

The Keratoconus table displays the following data:

- **KPI:** Keratoconus Probability Index

WARNING: Keratoconus scores are intended for use in a clinical exam and should never be used as the sole assessment tool in the diagnosis of Keratoconus.

The results of a clinical trial of the Keratoconus Prediction Index (KPI) are described in reference 1 and more information about mathematics are in reference 2.

- **SI:** symmetry index - the difference between the average thickness of two circular zones centered on the vertical axis;
- **AGC:** Apical curvature gradient - mean change per unit length of cornea power between the top of the cornea and the periphery of keratoconus;
- **AK:** Apical keratometry - curvature of the apex of the cornea.

Note: Values associated with a high probability of keratoconus are marked asterisk, like the SI and AGC values in this example:

Keratoconus	
KPI	37%
SI	*1.70
AGC	*1.68
AK	7.06

Eccentricity table

The eccentricity table shows the eccentricity of the cornea (P) in various predefined positions. The positions are determined by the meridian and the angle from the keratometric measurement axis. R0 is the apical position (apex of the cornea) and the angles are 10°, 15°, 20°, 25° and 30°. The meridians used are the nasal, temporal, lower and upper half meridians, as well as the horizontal and vertical meridians and the average of the horizontal and vertical meridians.

Eccentricity-P							
	Mer.	R0	10°	15°	20°	25°	30°
Nasal	354°	43.75 D	0.91	0.84	0.77	0.57	0.46
Temporal	174°	44.06 D	1.28	0.91	0.81	0.69	0.62
Inferior	264°	43.84 D	2.03	0.80	0.76	0.78	0.88
Superior	84°	45.92 D	0.56	0.43	0.41	-0.11	-0.11
Horizontal	174°	43.90 D	1.10	0.88	0.79	0.63	0.54
Vertical	84°	44.86 D	1.29	0.62	0.59	0.34	0.88
Medium	0°	44.38 D	1.20	0.75	0.69	0.49	0.67

Sagittal radius

The "Sagittal Radius" table shows the radius of curvature of the cornea at various predefined positions. The positions are determined by the meridian and the angle from the keratometric measurement axis. R0 is the apical position (apex of the cornea) and the angles are 10°, 15°, 20°, 25° and 30°. The meridians used are nasal, temporal, lower and upper half-meridians and horizontal and vertical meridians. Average value - the average value of the horizontal and vertical meridians.

Sagittal Radius							
	Mer.	R0	10°	15°	20°	25°	30°
Nasal	273°	7.79	7.89	7.90	8.07	8.28	8.52
Temporal	93°	7.88	7.90	7.97	8.06		
Inferior	183°	7.80	7.86	7.90	8.04	8.22	8.34
Superior	3°	7.77	7.86	7.90	8.05	8.13	
Horizontal	93°	7.84	7.90	7.94	8.07	8.28	8.52
Vertical	3°	7.79	7.86	7.90	8.05	8.18	8.34
Medium	0°	7.81	7.88	7.92	8.06	8.22	8.43

Meridian table

The Meridian table shows the radius of curvature of the eye in three steps: 3 mm, 5 mm and 7 mm, for the following values:

K1 & K2: two values of orthogonal radii, where K1 is the flattest meridian and K2 is

the steepest meridian

AVG: average value of K1 and K2 - average value of curvature between K1 and K2

CYL: cylinder measurement

Meridians						
	3 mm		5 mm		7 mm	
K1	7.90	96	7.91	97	7.95	93
K2	7.78	6	7.83	7	7.88	3
AVG	7.84		7.87		7.91	
Cyl	-0.66	96	-0.41	97	-0.40	93

SIM-K table

The SIM-K table shows corneal keratometry , for the following values:

K1 & K2: two values of orthogonal radii, where K1 is the flattest meridian and K2 is the steepest meridian

AVG: average value of K1 and K2 - average value of curvature between K1 and K2

CYL: Toroidal toricity - the difference in diopters between the curvatures of the principal meridians (K1 and K2) and the axis of the negative cylinder.

SIM-K		
K1	7.88	180
K2	7.88	90
AVG	7.88	
Cyl	0.00	180

Anterior chamber analysis


Subtab Anterior chamber The Analyzes tab of Cornea displays a chart and graph for anterior chamber analysis (pachymetry (for WAM 700+ model), corneal irido angle , corneal thickness). The user can edit this diagram by rearranging the measured points.

No.	Description
1	Editor Displays the editing screen. On this screen the measured points can be moved by the user.
2	Negative Displays the negative image of the cornea
3	Strengthen Improves image contours
4	Corneal thickness Displays corneal thickness values in multiple areas on a corneal image
5	Irido angles Irido values and positions Angle on the image.
6	ACD Displays the depth of the anterior chamber in millimeters.
7	W.T.W. Displays the distance from white to white (relative to the cornea).
8	Kappa corner Displays the angle between the center or cornea and the visual axis.
9	Cornea Scheimpflug image of the cornea. T Indicates that this is a temporary side eye. N Indicates that this is the nasal side of the eye. Orientation Indicates pachymetry orientation (for the WAM 700+ model) (in the example - a horizontal section). Central Corneal thickness (CCT) in microns is displayed in red.

Editor

The Edit screen allows the user to move the measured points in the image to get better values. Clicking on 1 of the three image sections will display a larger image. You can then drag each point to a new location.

Please note that if the measured results are not within the following limits, the edit screen will be automatically displayed when entering the anterior chamber analysis results for that eye

- $5^{\circ} < \text{Irido Angle} < 60^{\circ}$
- $1.5\text{mm} < \text{ACD} < 4.5\text{mm}$
- $8\text{mm} < \text{WTW}$
- Irido points must be in the appropriate EDIT window
- $-30^{\circ} < \text{Corner Angle} < +30^{\circ}$
- $200\mu < \text{Corneal thickness} < 800\mu$
- Corneal thickness values should be symmetrical about the center


Change the central thickness of the cornea

To edit the thickness of the central cornea:

Select the yellow cross and drag it to the outer border of the cornea.

Please note: this is only recommended if the recognition of corneal thickness is not satisfactory. This will change the values for the central corneal thickness and display only the CCT, other values along the periphery of the cornea will be erased.

Tonometry tab

The Tono tab displays the results of tonometric measurements, as well as the results of pachymetry. For each eye, the measurement results are displayed individually, and their average value is also displayed in the following columns:

essilor Smith John 07-Sep-2014 09:32:21

Measurement	IOP mm/Hg	Pachy	IOPc mm/Hg
1	11	539.1	11.6
2			
3			
AVG	11.0		11.6

Measurement	IOP mm/Hg	Pachy	IOPc mm/Hg
1	11	503.5	14.1
2			
3			
AVG	11.0		14.1

IOPC correction by: Ehlers 2

ACD	3.39 mm
ACV	212.9 mm³
Kappa α	7°
Decentration	0.4 mm

ACD	3.34 mm
ACV	181.0 mm³
Kappa α	3°
Decentration	0.2 mm

IC Angles	Left	Right
	45°	43°

IC Angles	Left	Right
	44°	47°


Summary Maps Cornea Tono Opacity Coeff. Aberr. Simulation CL Fitting Exit

No.	Description
1	Measurement Designation of measurement number and mean
2	IPO mmHg _ Displays the actual intraocular pressure that was measured.
3	Pachy Displays the thickness of the cornea at its center, in microns.
4	POc mmHg _ Displays corrected intraocular pressure values, taking pachymetry readings into account.
5	ACD Displays the depth of the anterior chamber
6	ACV Displays the volume of the anterior chamber.
7	Kappa corner Displays the kappa angle.
8	decentration Displays the decentration between the apex of the pupil and the visual point of the cornea
9	Angle IC Displays the angles of the inner camera for each eye.

Warning: The safety and efficacy of the device in measuring intraocular pressure in patients with more than 3 diopters of corneal astigmatism or more than 33 mmHg . intraocular pressure was not assessed .

Opacity


The Opacity tab graphically displays the results of several measures - ocular aberrations, corneal aberrations, retro-illumination, pachymetry (for the WAM 700+ model) - to give an overview of what the eyes of an opaque patient look like.

No.	Description
1	Ocular aberrations Displays values for high order ocular aberrations
2	Aromatization of the cornea Displays high-order corneal aberration values
3	Lens Displays lens retro illumination image and opacity
4	LOCS This subhead contains a menu of retro illumination images. After this tab is displayed, select the image that most closely resembles the patient's retroillumination image. See below for more details.
5	Transverse section Cross section of the cornea and inner chamber
6	Shack - Hartmann Hartman image

LOCS Screen

After selecting this tab, the following screen will appear:


No.	Description
1	Rear Select the image that most closely resembles the patient's retroillumination image. The image on top represents the softest back opacity, while the image on the bottom represents the heaviest back opacity.
2	Cortical Select the image that most closely resembles the patient's retroillumination image. The image on top represents the softest bark opacity, while the image on the bottom represents the strongest bark opacity.
3	Opalescence Select the image that most closely resembles the patient's retroillumination image 1 represents the least opacity, 5 represents the most opacity.
4	Color Select the image that most closely resembles the patient's retroillumination image. #1 represents the lightest color, #5 represents the darkest color.

Select the image that best matches the patient's cornea.

There are 4 categories to choose from: Back, Cork, Opalescence and Color.

Please note that only 1 of these 4 options can be selected.

After exiting the results screen, you will see a message:

Click "Yes" to save or "No" to discard the changes. The LOC selection can only be saved after retrieving the results, not immediately after the measurement.

Aberration coefficients

Aberration tab Coefficients » Zernike coefficient values are shown in units of mm for:

- eyepiece;
- cornea;
- internal.

There are three options for displaying aberration coefficients:


Main values: shows a selection of the most common aberrations; defocus , astigmatism, general lower order (LOA), general higher order (TOA), coma, spherical aberration, 4th order astigmatism, trefoil, and tetraphile . Aberrations are shown in the form of values and axes.

LOA: Shows the values for the low order aberration, total LOA, tilt, defocus , and astigmatism.

HOA: Shows values of higher order common aberrations, common HOA, trefoil, coma, tetratrium , 4th order astigmatism, spherical aberration, pentafoil , fifth order trefoil, and fifth order coma.

Eyepiece


Ocular tab shows the value of the Zernike coefficients of the whole eye from the wavefront measurement. Can display day/night data or right/left data.

No.	Description
1	Aberration coefficients Select displayed aberration coefficients
2	Left / right display Select right/left display
3	Photographic film / Mesopic display Select day / night
4	Graphical representation
5	Numerical representation
6	Scale Changes the scale of the graphical representation

Cornea


Corneal tab shows the values of the Zernike coefficients from the anterior corneal surface from corneal topography measurements. Can display day/night data or right/left data.

No.	Description
1	Aberration coefficients Select displayed aberration coefficients
2	Left / right display Select right/left display
3	Photographic film / Mesopic display Select day / night
4	Graphical representation
5	Numerical representation
6	Scale Changes the scale of the graphical representation

Interior

Internal tab Coefficient of Aberrations presents the differences between corneal aberrations and ocular aberrations for each eye. Grain coefficients are displayed for Corneal, Internal and total aberrations. For example, the graph below shows that most of the HOA, Defocus, and RMS are due to internal eye aberrations.

No.	Description
1	Photographic film / Mesopic display Select day / night
2	Aberration coefficients Select displayed aberration coefficients
3	Left / right display Select right/left display
4	Graphical representation
5	Numerical representation
6	Scale Changes the scale of the graphical representation

Modeling


Visual acuity measures the ability of the eye to discern object details and shape using HOA (High Order Aberration) and LOA (Low Order Abbreviations).

Visual acuity measurement applicable to topography (cornea) and wavefront (Eyepiece). The results for both options are interpreted in the same way.

Modifications available:

eyepiece
cornea


Eyepiece

No.	Description
1	Photopic selector Selects daytime display
2	Mesopic Selects night display
3	VA/PSF selector Allows you to switch between visual acuity map display and PSF (Point Spread Function) Display
4	MTF selector Displays the MTF chart.
5	Correction selector Select HOA to access the same image with HOA correction. Select LOA to get the same image with LOA correction. You can select Correction or Full Correction.
6	R/L Displays two cards, 1 for each eye. Both maps are displayed with the same option display.
7	NOA / LOA Displays two maps for each eye, 1 with HOA and 1 with LOA correction.
8	Day Night Image correction for day and night conditions for 1 eye.
9	Visual acuity of the image Reflects the actual quality of the patient's vision.
10	PSF Image of the cornea PSF (Point Scattering Function).

Cornea

Subtab The Corneal is identical to the Ocular as the measurement of visual acuity in relation to topography and wave front. The results for both options are interpreted in the same way.

No.	Description
1	Photopic selector Select daytime display
2	Mesopic Select night display
3	VA / PSF selector Allows you to switch between visualization map display and PSF display
4	MTF selector Displays the MTF chart.
5	Correction selector Select HOA to access the same image with HOA correction. Select LOA to access the same image with LOA correction. You can select Full correction or no correction at all.
6	RL Display two cards, 1 for each eye. Both maps are displayed with the same display options
7	HOA LOA Displays two maps for each eye, 1 with HOA and 1 with LOA correction
8	Day Night Image correction for day and night conditions for 1 eye
9	Visual acuity of the image Reflects the actual quality of the patient's vision.
10	PSF Image of the cornea PSF (dot spreading function)

CL installation


On this tab, you can place contact lenses on the patient according to the test results. The drop-down menu of the CL Fitting tab has several types of contact lenses listed according to their manufacturer, as

well as their sizes. Select a specific lens type and the software will calculate an approximate fitting for you.

Note:

The contact lens sizes listed in the software may have changed since this version was updated. To avoid miscalculations, we recommend checking all sizes directly with the contact lens manufacturer before prescribing contact lenses.


Lens

No.	Description
1	TLT (Tear Layer Thickness) Allows you to set the value by which the lens is raised in the center in relation to the cornea. When the TLT is changed, a new model lens or a new manufacturer must be selected; lenses that meet the new application criteria are displayed. The program "remembers" the change in the TLT setting and uses this criterion the next time you access the contact lens environment.
2	Correction Angle The K1 orientation allows you to choose the direction of the flattest meridian of the toric lens. The default value is the direction of the flattest Sim K meridian. You can change this setting to simulate lens rotation. This parameter must be set before selecting a manufacturer.
3	Manufacturer Choose a contact lens manufacturer
4	Design Select the most appropriate lens from the selected selection.
5	Diameter Select diameter lens
6	base curve Select Lens Base Curve
7	Display Displays an image of a contact lens superimposed on an image of the eye

8	Scale The scale used for the display image
---	--


Values

No.	Description
1	Partition display Displays the cross section of the lens
2	Displaying values


	Choose to display values on the contact lens map
3	Table of values Displays the values of the point at which the user clicks on the image
4	Refractive indices Displays refractive values for the patient
5	Keratometry values Displays patient keratometry values
6	Display Displays an image of a contact lens superimposed on an image of the eye
7	Scale The scale used for the display image

Display

No.	Description
1	Selecting the first scale
2	Choice of the second scale
3	Choice of the third scale
4	Pupil Displays the pupil on the map
5	Kerato Displays keratometry values and a cross on the map
6	corners Displays corners on the map
7	Refractive indices Displays refractive values for the patient
8	Keratometry values Displays patient keratometry values
9	Display Displays an image of a contact lens superimposed on an image of the eye
10	Scale The scale used for the display image

Actions

No.	Description
1	Move lens Moves the lens to the image of the eye
2	oblique lens Tilts the lens
3	Reset Lens Revert all changes
4	Refractive indices Displays refractive values for the patient
5	Keratometry values Displays patient keratometry values
6	Display Displays an image of a contact lens superimposed on an image of the eye
7	Scale The scale used for the display image

Patient Management

You can manage the patient database - add patients, modify patient information, search for patients, and delete patient records - on the Patient screen.

Adding new patients

You can add a new patient to the database before or after exams. You can also add multiple related patients to the database quickly by duplicating the current patient record and changing the parts that are different.

To add a new patient before taking an exam:

On the main screen, select the "Add Patient" option

- OR -

On the Patient screen, select the Add Patient button.

The Patient Information dialog box is displayed over the Patient screen

Add Patient Dialog Box Over Patient Screen

- Fill in the fields.
- Choose from the following options:

No.	Description
1	The note Select Note to add a note to the patient file. This note can be accessed by clicking on the symbol on the result screen, the measurement screen, or the patient screen. After adding the information, the note will look like this.
2	Cancel Cancels saving a new patient.
3	Save Save the new patient to the database, close the dialog box and display the patient screen.

4	Save and open Save the new patient to the database and leave the dialog box open. You can add an additional patient by editing the fields as needed and saving additional patient information. Note. This is useful if you want to create several new patient records at once. For example, you can use this option to create new patient records for all family members.
5	Save and edit Save the new patient to the database, close the Patient screen, and open the Measurement screen. This allows you to go directly from adding a new patient to the database to taking measurements on them.

Adding a patient after an exam

You can complete the exam without first providing patient information. When an exam is completed, the exam results are stored in the database as temporary records under automatically generated names - with "Temp" as the last name and a number as the first name. Temporary entries remain in the database until you manually remove them or remove the system during routine maintenance. At any time before deleting them, you can convert them to permanent patient records by adding patient information to them, or shortening the results of a temporary patient record and inserting them into a permanent record (see Moving Test Results to Another Patient).

Adding related patients

To add patients from the same family to the database:

- In the Patient Information dialog box, complete the fields for the first family member.
- Select Save and New. The patient is saved in the database, but their information is still displayed in the dialog box.
- Change the information in the dialog box. Select Save and New if you want to add additional family members to the database or either Save or Save and Diag if there are no additional family members. For more information about these options, see Adding a Patient Before an Exam.

Adding a Patient from an XML Database

The machine supports XML, allowing you to import data from a compatible XML database. To import data to the device, click the icon located in the upper left corner of the measurement screen. After importing the data, continue the measurement. The dimension will be saved with the name and data from the database.

If the icon is not available, it means there is no data to import. Start importing data, make sure your settings are set correctly in Export configuration data.

Changing Patient Information

To change patient information:

- On the Patient screen, select a patient from the patient list.
- In the patient list, select Edit Patient. The Patient Information dialog box opens.
- Change the fields if necessary.
- Select Save. The changes are saved and the dialog box closes.

Patient Selection


There are two ways to select patients: individually or as part of a group. Patients must be selected individually to view or change their personal information, and to view and manage their test results. Patients can be selected as a group to remove their records from the database.

To select a single patient to see patient information and test results:

On the Patient screen, select a patient from the patient list. The patient entry in the list is highlighted with a blue frame.

To select 1 or more patients to delete them:

On the Patient screen, in the patient list, check the box next to each patient you want to delete.


No.	Description
1	Checkbox selected (to select patients in a group)
2	Selected patient (as an individual)

Patient search

You can search for a patient by filtering the display to only show patients whose first and/or last names begin with the specified letter sequence.

To filter the list of patients:

- On the Patient screen at the top of the screen, in the Last Name field, enter 1 or more letters from the beginning of the last name you want to see. The list displays only those entries whose last names begin with the sequence of entered characters.
- In the "Name" field, enter 1 or more letters from the beginning of the name you want to see. The list displays only those entries whose names begin with the sequence of characters entered.

No.	Description
1	Filtered List
2	Last name filter
3	Name filter

Deleting Patients

To remove patients from the patient list:

- On the Patient screen, check the box next to each patient record you want to Delete.
- At the bottom of the screen, select Delete Patient. You will be asked to confirm that you want to delete the selected patients.
- Select Yes. Patient records are removed from the list.

Carrying out diagnostics

You can start a new diagnostic procedure by opening the Measurement screen. The Measurement screen can be opened from any of the following screens:

- Home screen
- Patient Screen
- Results screen

After opening the Measurement screen, you can prepare the patient and select test options as described in the Preparing to Run Diagnostics section, and then select Go to start diagnostic tests.

Initiating diagnostics from the main screen

To start a new diagnostic procedure from the main screen, you only need to select Change screen option (#2 in the diagram above).

Please note that if you open the "Measurement" screen directly from the main screen, no patient will be selected. The system automatically generates the patient's name for the results, with "Temp" as the last name and a number as the first name (for example, "Temp 1184"). Temporary entries remain in the database until you manually remove them or remove the system during routine maintenance. At any time before deleting them, you can convert them to permanent patient records by adding patient information to them (see "Changing Patient Information") or by shortening the results of a temporary patient record and pasting them into a permanent record (see Moving Test Results to Another patient).

Initiating diagnostics from the patient screen

From the Patient screen, you can initiate a new diagnostic procedure in several ways:

- Select "New Diagnosis" without selecting a patient record. The measurement screen opens with the patient's temporary name ("Temp" + number, see above).
- Select a patient record and select New Diagnosis. The Measurement screen opens for the selected patient. After starting the diagnostic procedure, the results are automatically saved in the record of the selected patient.
- For a new patient, in the Add Patient dialog box, fill in the patient information and then select Save and Diag. The patient record is added to the database and the measurement screen is opened for the new patient. After starting the diagnostic procedure, the results are automatically saved in the new patient record.

Initiating Diagnostics from the Results Screen

When the result screen is open, either immediately after completing a diagnostic procedure or when results have been loaded from a patient record (see Loading Test Results), select Measure to open the Measure screen for the patient. You can then run some or all of the tests again, or run additional tests.

Preparing to run diagnostics

After you open the Measurement screen, there are a few things you need to do before you start diagnosing:

- Clean the headrest.
- Place a new sheet of protective paper over the chin or clean the chin.
- On the "Measurement" screen, in the list of available measures (No. 6 in the diagram above), select the diagnosis you want to perform. The tests required to perform the selected diagnostic are highlighted in the list of tests (#18 in the diagram above).

- In the upper right corner (#7 in the diagram above), select whether you want to run the diagnosis in the right eye, the left eye, or both.
- In the "Position" section, select "Adult" if the patient is an adult or "Child" if the patient is a child. The head and chin head of the device moves to the default position for the selected patient type.
- If you are performing a WF test and/or Tonometer test (for WAM 700+ model) and want to change the parameters of these tests (number of measures taken), on the "Measurement" screen, select "Parameters". In the dialog box that opens, change the settings as needed and then click OK to save the changes.
- Ask the patient to sit down, place their chin on the chinrest and lean their forehead against the headrest.
- If you need to fine-tune the position of the chin, select the up and down arrows (#12 in the diagram above) as needed to move the chin up or down.
- If you need to adjust the position of the device head, select the up, down, right, and left arrows on the lens's viewing edge (#3, #16, #17, and #10 in the diagram above) as needed to move the head up, down, right, or to the left.
- Ask the patient to look at the target (balloon). You should see his eyes in the lens image on the Measurement screen.
- Align the patient's pupil with the center of the target.

After completing all these steps, select "Go" to start the diagnosis.


6.2 Management of the diagnostic procedure

The diagnostic procedure is controlled on the "Measurement" screen. From this screen, you can select the diagnosis you want to perform, position the block head and chin so that the patient's eye is in the center of the camera, change the default test settings, start and monitor the diagnostic procedure.

- Centering the patient's eye
- Skip measurement during diagnosis
- Speeding up the measurement process
- Stop diagnostics before it is completed
- Performing tonometry

Centering the patient's eye


Before starting the measurement, make sure that the patient's eye is in the center of the chamber:

No.	Description
1	Camera view Shows an image of the patient's eye that appears in front of the camera.
2	direction arrows Click on the arrows to center the patient's eye in the camera window.
3	Cross camera The patient's eye must be in the center of the camera's cross.
4	Chin arrows Use the arrows to raise or lower your chin

On the Measurement screen, you can see the camera view (#1), make sure the patient's eye is in the center of the cross (#2). To center the eye, the arrows on the sides of the camera's view (#3) move up, down, left, or right. If you cannot see the patient's eyes in the camera window, you may also need to raise or lower the chin rest by pressing the arrows located on the right side of the measurement screen (#4).

If the eye is visible in camera view but not quite centered:

Click on the center of the pupil (from the camera - a white arrow), the device will automatically start around the pupil.

Note. To ensure maximum accuracy, if the examiner accidentally moves the head of the device so that the patient's eye is not centered in front of the camera, the device will move so that it is re-centered in front of the eyes.

Missing measurement during diagnostics

During diagnostics, the Skip button is located at the bottom right of the screen (#15 in the measurement screen overview above). Select the button to skip the current test and proceed to the next test if there are more tests left, or complete the diagnostic and display the results on the results screen.

If a test is skipped during diagnostics, the circle denoting it is marked with a red X.

No.	Description
1	Test skipped
2	Next test continues
3	Skip

Accelerating the measurement process upwards

After running diagnostics, you can select a measurement to speed up the process if necessary.

To speed up the wavefront measurement process:

On the "Measurement" screen, select "Quick Mode". Note: This only works for wave front measurement.

Stopping Diagnostics Before Completing Diagnostics

When a diagnostic is running, you can stop it at any time.

To stop diagnostics before it is complete:


On the "Measurement" screen, select "Stop"


Glaucoma assessment (for WAM 700+ model)

The device provides an assessment of glaucoma by measuring both intraocular pressure (IOP) using aerometric tonometry, and measuring pachymetry to assess the depth of the chamber of the eye. You can

run the Glaucoma assessment as a standall measurement or select ALL measurements where it is automatically enabled.

After the device completes the pachymetry measurement, it will descend to the tonometer mode. It uses the centering and focus obtained from the pasha measurement to center and focus the tonometer measurement. The following screen will appear.


1	A plus Select the "+" button to move the focus closer to the eye
2	Minus Select the "-" button to move the focus away from the eye
3	Focus of the horn cross  Focused, the cross looks like this.
4	force Fire When you are satisfied with the focus, you can select Force Fire to take measurement or wait for the device to take a measurement.


Ask the patient to focus on the black dot. To automatically take a measurement, the apparatus must be centered and focused. The cross should be focused and transparent behind the big black dot. By default, the blood pressure monitor usually takes the measurement automatically, but sometimes you may need to help it focus and focus.

How to focus the tonometer


1. If 1 circle is around the other, click on the gray circle to align both circles. Once in the center, the tone is automatically triggered.



2. If the gray circle is blurred. Press the "+" button until you focus, and then the tonometer automatically fires.

Troubleshooting

If your camera view looks like this	What to do
	If, the cross is not focused, or the cross is not visible, but only 1 white circle is visible. Press + continuously until the cross becomes focused. If the cross is blurred again, you will pass the focus, press the - button until it becomes transparent. If in focus, the t1 is automatically triggered.

	The circles are not centered, the patient was moving. Click the house icon to return to the starting point, and then continue to focus and center the cross and circles.
	The cross does not focus, press continuously until the cross becomes more focused. If the cross blurs again, you passed the focus, press the button until you are satisfied with the focus. Tono automatically fires when the focus is correct.

Focused and focused, but no dimensions?


If the circles and cross are focused and cleared but the device has not taken a measurement, you can force the device to take a measurement by pressing Force fire . If the device stays in 1 place and doesn't seem to autofocus, check Tono Automation to see if the device is set to Semi .

Semi-mode

By default, the device has auto focus, which is fast and accurate. However, another focus mode is available in Semi mode . In this mode, the tonometer will remain at the focal point found with the pachymeter and focus correction must be performed by the examiner. Follow the steps above to set the focus. To select tonometer modes, go to Tono automation . **Note that in this mode the device will remain at the set focus point and all focusing will be performed by the examiner, no t1 measurements will be taken if the device is not in focus.**

Demo / Change number of measurements

If you would like to demo a tonal layer for your patient, select the Options button and then click Demo Puff . This pop-up menu also allows you to change the number of times each eye is measured by the tonometer.

Manual focus

Manual focus is automatically activated if the device cannot get good focus, allowing the examiner to manually correct the focus using the and buttons. It can be activated in the following modes:

- Topographies
- Pachymetry
- Tonometry

How do you know if you are in manual focus?

The tonometer always has manual focus, as well as autofocus, that is, its normal mode. However, other measurements have indicators if the device is not autofocus. Topo and Pachy measurements have an indicator:

Topo manual focus indicator :

During topo manual focus, the top left window on the measurement screen has this image:

Manual focus indicator pachymetry :

During the pachymetry focal point, you will see two indicators, the image above (with a home icon), as well as a warning message: "Step 1: Center and focus the rings and select Acquire."


Focused Tonometry (for WAM 700+)

Tonometry focusing is always available along with t1 autofocus. Instructions for

while improving t1 AF, move on to glaucoma shielding.

Topography focus


During Topography measurement, if the device cannot auto focus, it will switch to Manual focus. The following screen will appear:

1	Minus Select button to move focus further away from the eye
2	A plus Select button to move focus closer to the eye
3	Corneal rings These rings should be as sharp and as dark as possible.
4	Apply Press this button once you are satisfied with the focus of the cornea rings

Focusing Pachymetry

Just before the device switches to Pachy Manual Focus , the message "Step 1: Center and orient the rings and select ACQUIRE" will appear. Select "OK" and follow the steps listed in Topo Manual focus , Once this step is completed, you will see the message "Step 2: Use the focus buttons to align the Apex cornea , and then select Purchase." After selecting OK, you will see the following screen:


1	Corneal angle The top of the cornea should be centered across horizontally, vertically it should be slightly above the red cross (as seen in the image above).
2	Cross The cross is used to center the corneal apex vertically and horizontally.
3	Minus button Select the button to lower the top of the cornea
4	plus button Select button to move cornea elevation
5	Apply Completes the measurement

Sometimes the cornea is not centered over the cross (#2 above) exactly and may appear below or above it. For example, if it looks like this:

Close vision test

To add Near Vision Testing to Wavefront diagnostic :

On the Measurement screen, select the Options button:

1	Demonstration Choose to showcase
2	wave front Number of selected wave front parameters
3	Tono measurements Number of measured t1s
4	closest distance Changing distance from study Near Vision , can be adjusted from 30-60cm
5	Cancel Cancels the near vision test
6	Save Saves your preferences

6.3 Working with test results

The test results are displayed automatically on the results screen as soon as the diagnostic procedures are completed. You can choose from a variety of ways to display your results on the Results screen. In addition, you can print the results, export them to a file on an external device, or initiate additional diagnostics from the results screen. Test results are also automatically saved to the patient database so you can see them again at any time.

- Return / add dimension
 - Export results to Phoropter
 - Export results via network or USB
 - Export results to remote access
 - Print test results
- Saving test results

Return / add dimension

Remeasurement

If for some reason you are unhappy with the results of your examination, you have the option to cancel the measurement before exiting the results menu. Immediately after completing the measurement, select the "Measure" button. After selecting the diagnostic you wish to perform, the following screen will open:

If the selected

new diagnostic contains measurements that have already been taken, the device will ask you if you want to repeat them. Select the measurements you want repeated and then select "OK". Proceed with the measurement as usual.

Adding dimensions

To add measurements to the result screen, press the Measure button immediately after completing the measurement. If the measurements will not be repeated, this pop-up menu will not be opened. Select a diagnostic and continue with the test. Additional information will appear on your results screen along with your previous measurements.


Note. When you repeat a measurement using the "Measure" button, the previous results are erased.

Export results to Phoropter

On the "Results" screen, you can export the test results to the phoropter, which is either connected directly to the device via its serial port or USB port, or can be accessed over the network to which the unit is connected.

To export test results to a phoropter :

- On the "Results" screen, select the "Export" button.
- In the "Data transfer" section, select "External device".
- Select Submit. Test results are sent to the phoropter


Note. Make sure a template is selected in the Export Settings menu. it is properly connected to the machine.

Export results to a computer or USB stick

On the Results screen, you can export the test results to a computer that can be accessed. using the machine via LAN, or to a USB media (such as an external hard drive, flash drive, or memory card in a card reader connected to the USB port device) connected to the machine from the Export dialog box.

To export test results to network or USB:

On the Results screen, select the Export button. The Export dialog box will open

Export dialog

- In the "Data transfer" section, you can select your preferences, "Network" or "USB".
- Select Submit and this browser will open:


- Find your preferred location and select OK.

Note. Make sure the devices are connected to the machine

Export results to an external device

To export test results to an external device:

On the Results screen, select the Export button. The Export dialog box will open.

- In the "Data transfer" section, select "External device".
- Select Submit.

Note: Make sure that the devices are connected to the machine and that the external device is set up correctly in the export settings.

Export results to email

On the Results screen, you can export the test results to the default recipient's email or to a new recipient.

To export test results to an email message:

On the Results screen, select the Export button. The Export dialog box will open.

Data Transfer

☐ External Device

☐ Network

☐ USB

☐ Email

Print

☐ Ticket

☐ Print Screen

Report

☐ Print

☐ Preview

Cancel Send

Select Email and SEND and the following screen will appear:

The following screen will appear:

From 1: [www.goszdravnadzor.ru]

To 2: [meira, Weird, Colleague 3, Patient, Laboratory 1, lets, Laboratory 3]

Subject 4: TEST RESULTS

Message 10: THIS 5 CONFIDENTIAL em 6

Set as default

Set as default

Import Note 11

Cancel 13 Send 14

No.	Description
1	Where By default, the sender's name is displayed from the settings menu.
2	Where Enter the email address you would like to send the information to in this field, or go to page 3 and select a default recipient.
3	Default Recipients

	The default recipients are displayed here. To send an email, select the checkboxes next to the names you want to email. Both Peers and Laboratories are available for selection. To set up default recipients, go to Print / E-mail.
4	Subject It can be free text or a default object. To select the current default object, press the 9, Set Default button.
5	Summary Select to include a summary screen in the message.
6	Topography kit Includes h8maps.
7	Patient File This is the patient data in a zip file.
8	Report This is a patient report pdf file.
9	Set as default Sets the theme as the default object
10	message box Add your message here.
11	Import entry Includes a note from the patient file in the body of the message.
12	Set as default Sets the message as the default message.
13	Cancel Cancels the message.
14	Send Click to send a message.

Note: In order to send an email, you first need to configure the email settings in the "Print / Email" section of the settings menu.

Printing test results

From the Results screen, you can print test results by exporting them to the built-in printer.

To print test results using the built-in printer:

- On the "Results" screen, select the "Export" button.
- In the "Print" section, select "Ticket".
- Select Submit.

The dialog box contains three sections with checkboxes:

- Data Transfer**
 - ☐ External Device
 - ☐ Network
 - ☐ USB
 - ☐ Email
- Print**
 - ☐ Ticket
 - ☐ Print Screen
- Report**
 - ☐ Print
 - ☐ Preview

At the bottom are two buttons: **Cancel** (with a close icon) and **Send** (with a folder icon).

To print a screen capture of your results:

- In the "Print" section, select "Print"
- Select Submit

Note. To print to an external printer, you must select it in the menu's print export settings.

Create a report

To create a report:

- On the "Results" screen, select the "Export" button.
- The following dialog box will open:

This is a duplicate of the dialog box shown above. It contains three sections with checkboxes:


- Data Transfer**
 - ☐ External Device
 - ☐ Network
 - ☐ USB
 - ☐ Email
- Print**
 - ☐ Ticket
 - ☐ Print Screen
- Report**
 - ☐ Print
 - ☐ Preview

At the bottom are two buttons: **Cancel** (with a close icon) and **Send** (with a folder icon).

In the "Print" section, select "Report".

- Select Submit.

The report viewer window will open:

No.	Description
1	Template menu Choose a report template according to your preferences.
2	It shows a preview of the maps and data available in the report
3	Select "Print" to print the report, or "Export" to send it to another device.

For more information about the possible report configurations, see the "Report Setup" section

Saving test results

The test results are automatically stored in the database. If a patient was selected before starting the diagnosis, the test results are stored under the patient's name. If no patient was selected at the start of the diagnosis, the system generates a temporary name with the last name "Temp" and a number for the first name (for example, "Temp 1185").

Temporary test results are automatically purged from the system at set intervals (once a day, week or month, depending on configuration settings, see "Purge History"). If you want to keep a temporary test result permanently, you must either cut the results from the temporary record, paste them into a normal patient record, or add the correct patient information to the temporary record. Both of these actions can be performed from the Patient screen (see Moving Test Results to Another Patient and Changing Patient Information).

Managing the list of test results

Stored test results are managed on the Patient screen. The results are stored in the name of the patient.

To view a list of test results for a patient:

- On the Patient screen, on the left side of the screen, select a patient record. the patient's record is highlighted with a blue box and their test results are listed on the right side of the screen.

- Select test results
- Download test results
- Export test results from the patient screen
- Transferring test results to another patient
- Deleting test results

Selection of test results


There are two ways to display test results that are stored in the database: individually or as part of a group. To select them individually, select the result you want to change and then you can view, print or export it. Results can be selected as a group to be moved from 1 patient record to another or removed from the database.

To select 1 test result:

- On the "Patient" screen, select a result from the list of results. The result is highlighted with a blue frame.

To select 1 or more test results to move or delete:

- On the Patient screen, in the list of results, check the box next to each result. want to delete.

No.	Description
1	Featured Patient
2	Check box selected (to select results as a group)

3	Result selected (as separate)
---	-------------------------------

Download test results

To download test results:


- On the "Patient" screen, select a result from the list of results. The result is highlighted with a blue frame.
- At the bottom of the list, select Download. The selected test result opens in the results screen.

Export test results from the patient screen

Test results can be exported from the patient screen to the built-in printer for printing as a ticket or as a file on an external device connected to the device.

To export test results from the Patient screen:

- In the "Measures" list (#6 on the "Patient" screen above), select the measure, export. The measure is highlighted in blue.
- Select "Export" at the bottom of the list. The Export dialog box will open.

- To print the test results on the built-in printer, in the "Print" section, select "Ticket" and then "Submit". The ticket is printed.

- OR -


To export test results to an external file, in the Transfer section, select an external device and click Submit. A dialog box will open with a file splash screen. Navigate to the folder where you want to save the file, enter a name for the file, and select Save.

Export test results from the results screen

Test results can be exported from the results screen to the built-in printer for printing as a ticket or as a file on an external device connected to the device.

To export test results from the results screen:

- Select the "Export" button. The Export dialog box will open.

No.	Description
1	External device Select to transfer data from the device to an external device; for example, XML data (and XML databases) or data into a phoropter head . View external device, go to the export settings menu.
2	Net Select to transfer data from the device to the network. To view the data options and configure the export path, go to the export options setup menu.
3	USB Select to transfer data from the device to USB. To view the data export options, navigate to the export settings menu. Warning - before connecting the USB to the device, make sure that it has been screened for viruses for the first time.
4	Email address Select to transfer data to email recipients. Go to Email Settings to set up the host and recipient options.
5	Ticket Allows you to print ticket data on the device's internal printer
6	Screenshot Select to print a screenshot of the summary screen, this requires an external printer connection.
7	print Prints the report selected by default in the Print menu, this requires a connection to an external printer.
8	Preview Opens a preview of the measurement results report
9	Send Confirm selection and transfer

- To print the test results on the built-in printer, in the "Print" section, select "Ticket" and then "Submit". The ticket is printed.


OR -

To export test results to an external file, in the Transfer section, select an external device and click Submit. A dialog box will open with a file splash screen. Navigate to the folder where you want to save the file, enter a name for the file, and select Save.

7.5 Transferring test results to another patient

To transfer test results to another patient:

- On the "Patient" screen, select a result from the list of results. The result is highlighted with a blue border.
- At the bottom of the list, select Cut. The selected test result is removed from the current patient record. The patient selection is cleared and as a result no patient is selected and test results are not listed in the results list.
- Select the patient record to which you want to move the test results. The patient record has a blue border and the Cut button has been replaced with a Paste button.

- Select Paste (No. 1). You will be prompted to confirm that you want to insert the selected test results into this entry.
- Select Yes. The test results are added to the selected patient record and displayed in the test results list on the right side of the screen.

7.6 Deleting test results

To remove test results from a patient record:

- On the Patient screen, check the box next to each test result you want to delete.
- At the bottom of the screen, select "Delete measure". You will be prompted to confirm that you want to remove the selected measures.
- Select Yes. The test results are removed from the patient record.

Database Management

The patient database contains all patient information stored in the system. All test results are linked to specific patient records; even if you have previously taken diagnostic tests without entering any patient


information, the results are stored in the patient database in automatically generated temporary patient reports. You can also export the entire contents of the database to an external file as a backup; import patient records that were previously exported; and delete all patient records or all temporary patient records from the machine database.

Database export and import


You can export the entire database to a folder stored on an external device such as a USB hard drive or flash drive, or on a PC that is on the same LAN as the machine. You can also import data from any manufacturer's database. You can either import the entire database or just the data for 1 patient.

To export your database to an external device:

When you select the Maintenance tab, the following screen appears:

in the "Database" section on the right, select the "Backup" button



The Browse Folder dialog box opens.

- Navigate to the device where you want to save the exported database data and select (or create and select) the folder where you want to save the data.
- Select OK. The backup starts. A ClientDB folder is created in the selected folder and the data from the database is stored in it. The "Copy" dialog box will open and show the progress of the backup, and the message "Backing up the database" will appear.

When the backup is complete, a new message will appear asking if you want to clear the data from the machine's database.

- If you want database data machine,
Otherwise, select No.

to delete all
on the
select Yes.


Note: For more information about deleting all data from the database, see Deleting all patient records.

To import data from an exported database:

- On the "Configuration" screen in the "Main10ance" tab in the "Database" section select Import Data. The Browse Folder dialog box opens.
- Navigate to the device where the database data resides and select the ClientDB folder where it is stored.

Note: If you wish, you can import data for a single patient rather than all the data in the database. To do this, in the ClientDB folder, select the patient folder. (The folder name contains the patient's full name, date of birth, and gender. For example, a folder called Melodie # Emma # 12-03-65 # f contains data for Emma Melody , whose birthday is 3/12/65 and whose gender is female.)

- Select OK. The data import process begins. The Copy dialog box opens. import progress, and the message "Restore database in progress" is displayed. When the import process is complete, the dialog box and message will disappear.

Note: During the import process, you may be asked to confirm that you want to overwrite the existing folder with the imported 1 if both have the same name. Data for new patients will be added to the existing database and will not be overwrit10 in your database. Select "Yes" to overwrite existing data.

Deleting all temporary patient records

Temporary patient data is automatically deleted from the database every day, week, or month, depending on the configuration settings (see Device Settings). In addition, you can manually delete all temporary data at any time.

To remove all temporary patient records from the database:


- On the Patient screen, in the Last Name filter, enter Temp . Only patient records last name fields start with Temp . This includes all temporary patient records in the database.
- In the title bar of the patient list, check the box (#8 on the patient chart screen).
- Ensure that there are no patient records that are not temporary records (i.e., records of patients whose last names begin with " Temp "). If yes, please clear the receipt boxes next to their names.
- In the Patient list, select Delete Patient. You will be asked to confirm that you want to delete the selected patient records.
- Select Yes. The selected patient records are deleted.

Deleting all patient records

If necessary, you can delete the entire patient database.

To delete the entire patient database:

- On the "Configuration" screen in the "Main10ance" tab in the "Database" section select Delete Data. You will be asked if you want to clear the data from the machine.

Select Yes. You will be asked to confirm that you want to delete all entries in the database.

- The Uninstall dialog box opens and shows the progress of the uninstall process. When the process is completed, the dialog disappears.

General actions and functions

This chapter provides more information about the machine and how to use it: using the screen saver and password-protecting the screen, entering text in text fields, connecting a keyboard and mouse to the machine, updating the software, and getting help.

Launching the screen saver

You can launch the screen saver from the Home screen or the Patient screen.

To start the screen saver:

- select the "Screensaver" button.

The device can also be configured to launch a screen saver automatically after a set number of minutes without any user interaction (see "Screen Saver").

Note: We recommend adding a password to the splash screen to prevent unauthorized users from accessing it. To learn more about adding a password, go to the Screen Protector section.


Screen protection

You can protect the system by requiring a password to disable the screen saver. Ideally, you should also set the screen saver to start automatically after a few minutes without user activity. See Screen saver for more information.

Entering text

Whenever you select a text field, such as the field where you enter a patient's name (see "Adding New Patients"), a virtual keyboard will open on the screen allowing you to enter characters. Also, if you wish,

You can connect a mouse and/or keyboard to the device's USB ports. If you connect a keyboard, you can enter characters into text fields using the virtual keyboard or the physical keyboard. Similarly, if you connect a mouse, you can select characters on the virtual keyboard by tapping them on the touch screen or by clicking the mouse.

To close the virtual keyboard:

Select a  button.

Software update


When the manufacturer releases a software update, you can use it to update the software on the unit.

To update the software:


- After receiving the update file from the manufacturer and saving it either on the machine itself or on a device to which it can connect. (For example, download the update file to a USB flash drive, and then connect the flash drive to the device's USB port.)

From the Home or Patient screen, open Configuration screen.

- Select the Main10ance button. The Main10ance tab opens.
- On the left side of the screen, under Software Main10ance, select Software Upgrade.

A file selection dialog will open

- Navigate to the update file.
- Select OK.

Integrated Help

You can access the machine's Help Guide directly from the software. Click the icon to access help files. To close the help files, click the "X" located in the upper left corner of the  screen.

Getting Help

If your device is connected to the Internet, technical support staff can connect to the device online as needed to provide guidance and troubleshoot problems. This feature uses H2O, a software package that allows remote control and PC-to-PC sharing. You can activate H2O from the Configuration screen under the Main10ance tab. For more information about this setting, contact technical support.

Device setup

Access to configuration settings

To open the configuration screen from the main screen:


Select the Configuration Screen button (#5 in the diagram above).

To open the Configuration screen from the Patient screen:

Select the "Configuration" button (#12 in the diagram above).

Overview of the configuration screen

The Configuration screen gives you access to the machine settings and system main10ance functions.

No.	Description
1	Tab buttons Select the button to display a tab containing configuration options for the category, or select the Main10ance tab to access system main10ance functions such as software updates and hardware calibration.
2	Sub tab buttons Some tabs are divided into sub-tabs. Select a button to display a sub-tab.
3	Function buttons Standard buttons available on all configuration screen tabs: - Disable: turn off the device (see Turning the device on and off) - Remove USB: Safely remove a hardware device connected to the device's USB port. - Factory: Restore factory settings. (Note: This does not delete data.) - Cancel: Cancel all changes made to the configuration settings. - Save: save all changes made to the configuration settings, close the Configuration Screen and display the main screen.

Changing configuration settings

You can change any of the configuration options that appear on the various tabs of the configuration screen.


o change configuration settings:

- Change the settings you want to change. You can change settings on more than 1 tab without saving the changes, as long as you don't close the configuration screen before saving them.
- When you have changed everything you want to change, select "Save". Changes , and the previous screen is displayed.

General settings

The General tab of the Configuration screen is divided into three sections:

- Regional options: date, time and language settings
- Device Options: Choose a default screen and schedule how of 10 temporary data is deleted from the database
- Security: screen saver settings.


No.	Description
1	Regional Options
2	Default Screen Select the screen you want to use as the default screen - the screen that appears when you first turn on the machine and when you select "Exit" to close the "Results" screen. The default screen can be the home screen, patient screen, or measurement screen.
3	Device Options
4	Safety
5	Screen saver Check the box to set the screen saver to automatically activate when there has been no user activity for the number of minutes selected in the "Please wait" section (below).
6	Please wait (min) Displays the number of minutes of inactivity that will cause the screen saver to activate. Select the "+" and "-" buttons to change the number.
7	Password protection Check the box to require a password to disable the screen saver.
8	Change password Select this button to change your password.

	In the dialog box that opens, in the "Old password" section, enter the current password (there is no password by default.) In the "New password" section and under password confirmation, enter a new password. Click OK to close the dialog box.
9	Automatic revert to default Allows the user to exit the results screen and automatically return to the default screen after the selected time. Select this button to change the time according to your preference (eg 2 min, 3 min, Never).
10	Clear data history Select how often you want the system to delete temporary patient records from the database. Note: Temporary patient records are measurements that have been taken without being associated with a patient record. They have the last name "Temp" and a number for the first name.
11	warning signal Choose whether you want to hear a beep when a warning message is displayed.
12	Formatting time Select the clock display type to use for the rest of the interface: - 12 o'clock: 12-hour clock, with AM and PM to indicate which half of the day is included - 24 hours: 24 hour
13	minutes Displays the current minutes per hour. Select "+" and "-" buttons to reset current minute.
14	Watch Displays the current hour. Select the "+" and "-" buttons to reset the current hour.
15	Language Select interface language
16	Format date Select the date display format in the interface: - dd / mm / yyyy: European standard, for example June 26, 2013, will be displayed as 26/06/2013 - mm / dd / yyyy: American standard, for example June 26, 2013, will be displayed as 06/26/2013
17	the date Select the current date.

Diagnostics


The Diagnostics tab of the Measurement tab allows you to select which diagnostics will be available in the list of diagnostic functions on the Measurement screen. You can also choose which of the diagnostic tests will be the default - the diagnostic that will be selected when you first open the measurement screen.

To customize the list of diagnostic programs:

- In the Display (No. 1) column, check the box for each diagnostic test that you want to be available on the Measurement screen.
- In the "Default" column (No. 2), select the default diagnostic - the diagnostic that will be displayed when you open the measurement screen.

Scroll down for more diagnostics


Diagnostics	Description
All	Includes WF, Topo, Tono, Keratometry, Pachy and Refraction
Glaucoma	Includes Pachy and Tono measurements
Refraction	Includes refraction measurements
wavefront	Includes Wavefront, Aberrations and Refraction measurements
topography	Includes topography measurements
Pachy	Includes Pachy measurements
Cataract	Includes Wavefront, Topography and Pachy measurements
Corneal app.	Includes Wavefront, Topography and Pachy measurements
C.L. Fitting	Includes Wavefront and Topo measurements
Pupillo	This measurement shows pupil size and reaction time according to various lighting conditions.
Custom	This measurement can be customized according to your preferences and may include all measurements in the device except pupillo.

To select a dimension, select it from the list above. Only measurements selected in the default menu appears in the measurement menu.

Default

The default Tab lining on the Measurement tab has two sections. In the first section, Default Device Settings, you can set the default sound and brightness settings, and choose whether the default headrest and chin posts are for adults or children. In the second section "Measurement Defaults", you can select the default value that should be d1 by default in order to perform Wave FrontandTonometry tests and whether to include near vision in wavefront measurements. You can change the position setting and the number of times tests are to be performed "on the fly" on the "Measurement" screen before you start a measurement.


No.	Description
1	Horn Select the volume of beeps that the machine emits during the measurement. Measurement default settings - WF measurements: displays the default number of times before the Wave Front test should be repeated. Select the + and - buttons to change this value. - Tono Measurements: displays the default value specified in the Tonometry test should be repeated. Select the + and - buttons to change this value.
2	Target luminosity Select the brightness level of the fixation target. Topo with Ringer Display the topography results with the Ringer option after taking a measurement and before displaying the results.
3	Job title Select whether the default position of the headrest and chin should be suitable for adults or children.
4	Reading Distance Set the reading distance for near measurement, from 30cm to 70cm.
5	Add Dimension with WF Select to add Near Vision to the Wave Front test. Select + and - to change the default distance for near vision testing. Note: This will add Near Vision testing to each wavefront test. T1 Measurements Allows you to change the number of Tono measurements.
6	WF measurements Select the number of wavefront measurements for each wavefront exam. Only measurements that include the wavefront.
7	Tono measurements Tono Measurements: displays the default value specified in the test Tonometry should be repeated. Select the + and - buttons to change this value. Only exams that include tono will be completed.
8	Tono Automation Two types of focus are available for tonometer, auto and semi. By default, t1 measurement is in auto focus mode, where the device can take measurements with little

	intervention from the examiner. In Semi mode, the tonometer remains at the centering point found in the pachymeter measurement and all focusing is d1 by the examiner.
9	Topo with Ringer By default, the call is editable. To remove the shown ringer display immediately after topo measurements, clear this check box.

illumination

We recommend using the default lighting settings, although they can be changed in this screen:


No.	Description
1	Light level Use the + and - buttons to adjust light levels
2	Waiting time Use the + and - buttons to adjust the waiting time between light levels.
3	Save Select to save changes.

Results Settings

- Refractive abrasives
- W.F.
- Machine User Manual
- Device setup
- Topography
- Tono

Refraction abrasion settings


The Refraction and Aberration tab underlay of the Results tab contains a display of default settings for the refraction and aberration data of the results screen.

No.	Description
	 The screenshot shows the Essilor software interface. Numbered callouts point to the following elements: 1: Terminology section (Standard, Absolute, Auto, Normalized, Tangential, Instantaneous) 2: Refraction section 3: Step section (0.25, 0.12, 0.01) 4: Cylinder section (+, -, Auto) 5: Vertex Distance section (Dists., 12mm, 13.5mm, 13.75mm, 15mm) 6: Aberrations section 7: Zernike Unit section (Micron, Diopter) 8: WF Presentation section (Day/Night, 3/5) 9: Aberrations Aperture section (3.00mm, 5.00mm)
1	Terminology You can set the terminology for the map and scale to Tangential or Instant (for map) and Standard or Absolute, Auto or Normalized (for scale)
2	Refraction Default settings for refraction and ocular wavefront
3	Step Select the measuring step.
4	Cylinder Select a negative or positive cylinder convention, or select "Auto" to set the cylinder sign to match the sphere's equivalence sign. Select the measurement step.
5	Distance to summit Select the default vertex distance. (You can change this value on the fly in the Results Screen.)
6	aberrations Default settings for cornea and ocular aberrometry
7	Zernike block Select the unit of measurement (micron or diopter) for the measurements shown.
8	WF presentation Select Day/Night to represent refractometry measurements as photopic and mesopic values, or 3/5 to represent refractometry measurements as values for 3mm and 5mm student apertures.
9	Deflection Hole

If you selected "Day/Night for WF Presentation" (see above), select "Day" to present the default measurements in terms of photopic pupil aperture, or "Night" to present them in terms of mesopic pupil aperture.

If you selected 3/5 for the WF presentation (above), choose whether to present the default measurement results as a 3 mm pupil aperture or a 5 mm pupil aperture.

WF


No.	Description
1	Default map Choose which cards to display on the Results screen by default
2	Color palette The colors for maximum and minimum curvature are always the same in the color palette. The spacing size determines the remaining colors. - ANSI has 21 to 25 steps from red to blue. In particular, for ANSI Z80.23, red, orange, and yellow are steeper areas of the cornea. Green is the "Intermediate" color. The light and blue 1s show the flattest areas. - ISO has 26 steps from red to blue - The custom palette has 26 steps from white to blue. The central value is yellow. In the settings, you can select this color palette and apply green as the central value. Note. In case any value is greater than the maximum interval or less than the minimum interval, these areas will correspond to the color assigned to the maximum interval or the lowest interval, respectively.
3	Step Select a standard rate of change for the color ramp: value range represented by each individual color. Select "+" to increase the step or "-" to decrease it. Note: You can manually change the step on the result screen as needed. This parameter specifies only the initial setting. Note: When ANSI is selected in the Map section, the pitch is automatically set to 1.0 and the selection cannot be changed.
4	Central value Select the default value for the color scale. Select "+" to increase the value, or "-" to decrease it.

	Note: You can manually change the center value on the Results screen as needed. This parameter specifies only the initial value. Note: If ANSI or ISO is selected as the map, the center value is automatically set to 44.00 and the selection cannot be changed.
5	Green in the middle Select this option to orient the palette so that the center value is displayed in green instead of the default yellow. Note: If ANSI or ISO is selected for the map, this option cannot be selected.
6	Show by default Selecting Options to Display on the Results Screen

Topography

The Topography tab of the Results tab has two sections. The top part contains general settings. The Palette Topography section contains options that determine how topographic maps will be color-coded.

The color palette on the left side of the screen changes dynamically as you change the settings in the Palette Topography section. In addition, when certain options are selected, other settings in this section are automatically set and their fields become unavailable. For example, if you select ANSI in the Map section, all other options in the section are set automatically and cannot be changed.


No.	Description
1	Default map Choose which cards to display on the Results screen by default
2	Unit ism Select the default units of measurement - mm or diopters. Note: If necessary, you can change the units of measure on the results screen. This setting determines only the initial units used. Note: If ANSI is selected in the Map field, the Diopter field is automatically selected in this case and the selection cannot be changed.
3	standard scale At this scale, colors correspond to values to ensure consistency when compared between maps. There are several possible standard scales: ANSI, ISO, Klyve/Wilson and Maguire/Waring. You can set up a custom standard scale by setting the center value and pitch. Then the maximum and minimum values are automatically determined. Choose a standard color coding convention:

	- ANSI: Use the ANSI color coding convention. All other settings in the Topography Palette section are set automatically as required by convention and cannot be changed. The average will be set to 44D and step up to 1D. - ISO: Use ISO color coding. Scale, center value and color the central value are set automatically according to the requirements of this agreement and cannot be changed. The average value is set to 44D and the step can be changed to 0.5, 1.0 or 1.5D. - Custom: allows you to customize the palette as you wish; all other settings available. Note: These options determine the default options that are used in topographic maps on the results screen. You can always change the units, pitch, and center value from the Maps tab.
4	Step Select a standard rate of change for the color ramp: value range represented by each individual color. Select "+" to increase the step or "-" to decrease it. Note: You can manually change the step on the result screen as needed. This parameter specifies only the initial setting. Note: When ANSI is selected in the Map section, the pitch is automatically set to 1.0 and the selection cannot be changed.
5	Central value Select the default value for the color scale. Select "+" to increase the value, or "-" to decrease it. Note: You can manually change the center value on the Results screen as needed. This parameter specifies only the initial value. Note: If ANSI or ISO is selected as the map, the center value is automatically set to 44.00 and the selection cannot be changed.
6	Color palette The colors for maximum and minimum curvature are always the same in the color palette. The spacing size determines the remaining colors. - ANSI has 21 to 25 steps from red to blue. In particular, for ANSI Z80.23, red, orange, and yellow are steeper areas of the cornea. Green is the "Intermediate" color. The light and blue 1s show the flattest areas. - ISO has 26 steps from red to blue. - The custom palette has 26 steps from white to blue. The central value is yellow. In the settings, you can select this color palette and apply green as the central value. Note: In case any value is greater than the maximum interval or less than the minimum interval, these areas will correspond to the color assigned to the maximum interval or the lowest interval, respectively.
7	Green in the middle Select this option to orient the palette so that the center value is displayed in green instead of the default yellow. Note: If ANSI or ISO is selected for the map, this option cannot be selected.
8	Show by default Selecting Options to Display on the Results Screen


Note: To ensure maximum accuracy, if the examiner accidentally moves the head of the device so that the patient's eye is not centered in front of the camera, the device will move so that it is re-centered in front of the eye. If the maximum deviation of the decentrator is not met, the device will not take measurements.

Tono

Here you can choose which correction to use for Tonometry calculation. Correction formulas are given below.

It is also possible to select a custom offset in which the user can insert values for the A and B coefficients.

No.	Description
1	Tonometer corrections Choose between different correction methods
2	Coefficient A Insert values for factor A
3	Factor B Insert values for factor B
4	Corrected IOP The custom correction is calculated according to this adjusted IOP formula

Export settings


The export options define options that determine how the results of the selected (or current) diagnostic can be exported to other devices or printed.

Note: These settings apply to the processing of individual test results. For information about exporting the entire patient database, see Exporting and Importing a Database.

- Data export
- Print / E - mail
- Automatic actions
- Remote access
- Reports

Data export

The Data Export tab of the Export tab contains options that control what happens when you export current test results to an external storage device or optical instrument.

No.	Description
1	USB Several file options are available for USB export according to your preference, fill in the checkboxes of the data you want to export.
2	Net Several file options are available for exporting over the network according to your preference, fill in the checkboxes of the data you want to export.
3	Export catalog Select the path to the directory you want to export.
4	External device Check this box if you want to export it to an external device. External devices communicate with the LANBOX and XML protocols to create XML files sent over the network. This allows you to transfer files to EMR and other equipment such as APH 400, phoropter and Remote Access. The test results will be transferred to the device in the selected format. Remote access can be accessed through the URL path specified in the

	XML data for each patient. Once configured, any computer on the intranet can view the patient file remotely.
5	XML export After checking the external device, select this field to specify the data transfer path for the XML export.
6	XML import Select to import XML data (patient name and information) from the specified data path. Double click the field to select the location of the imported data, click SAVE. Once the XML import is set up, the data will be available for import on the measurements screen. For more information, see Adding a Patient via XML.

File export options

The settings in the first section determine the purpose of the exported files, their format, and the patient information they should contain.

The first line defines the export options for devices connected to the machine's USB port. (When you export a file to a USB device, the Save As dialog box opens, allowing you to navigate to a folder and specify a file name for the file.)

The second line defines the settings for exporting to devices to which the machine connects via LAN. You can specify the path to the device and folder on the device where the files should be saved.

To set file export options:

- To enable export to a device connected to the USB port, select the first checkbox in the USB line (i.e. the checkbox immediately to the right USB).
- To enable export to a device that can be accessed via the local network, select the first checkbox in the Network row.

For each option you choose (USB and/or network), select the following options:

Option	Description
Screenshot (JPG)	Select this option to export the data as a JPG image file.
Patient File	Select this option to export data to the machine. This creates an XPS file containing test results and patient information (as specified in Naming Patient Information) in the format used by the machine. Test results and patient information can be imported from this file into another database.
Surname	Select this option to include the patient's last name in exported data.
Name	Select this option to include the patient's name in exported data.
An identification number	Select this option to include the patient ID in the exported data.


- If you select more than 1 option, all the options you have selected are activated. For example, if you select both Screenshot and Patient as the JPG image of the result, the Summary tab and the Patient file will be exported.

- If you have activated the "Network" option, in the "Export" section, enter the path to the folder where you want to store the exported data. To speed up the export to the network, we recommend mapping a network drive, for more information go to How to create a network map Drive a car.


How to create a network drive map

To map a network drive:

1. Using an external keyboard, press CTRL + F4
2. Select the application "Exit"
3. Press + e  to open Windows Explorer.

4. Select Map Network Drive

5. In the Drive field, select a letter to represent the new network drive
6. Select "Browse" and select the folder you want to access from the network
7. Select "OK" and then "Finish".
8. Now this driver is available for selection in the export pop-up window and in settings-> Export menu

Print / E-mail

The Print tab of the Export tab contains options that control what happens when you choose to export test results to a printer for printing. On the left side of the screen, ticket printing is configured - printouts using

the built-in printer of the device. The right section configures printing performed by an external printer connected to the device's RS232 serial port or 1 of its USB ports.

Note: To install an external printer, contact technical support or your agent.

The screenshot shows the 'Essilor' configuration window. On the left, the 'Ticket' section has checkboxes for 'Ticket' (checked), 'Ticket Header' (checked), and 'Print' (unchecked). Below these are input fields for 'Title text' (containing 'Dr. Smith Office') and 'Printer' (a dropdown menu). The 'Format' section has radio buttons for 'Summary', 'Normal', and 'Complete'. The 'Character Size' section has radio buttons for 'Large' and 'Small'. On the right, the 'Email' section has a 'Test Email' button and checkboxes for 'Email' (checked), 'From' (containing 'JohnSmith@company.com'), 'Host' (containing 'mail.company.com'), 'Password' (masked), and 'Port' (containing '16'). There are also checkboxes for 'SSL' and 'TLS'. Below the email section are sections for 'Colleague' and 'Laboratory', each with 'Name' and 'Address' fields. At the bottom, there is a 'Maintenance' section with a 'Turn Off' button, a 'Test USB' button, a 'Factory' button, a 'Cancel' button, and a 'Save' button. Numbered callouts 1 through 16 point to various elements: 1 points to the 'Ticket' checkbox, 2 to the 'Ticket Header' checkbox, 3 to the 'Title text' input field, 4 to the 'Format' radio buttons, 5 to the 'Character Size' radio buttons, 6 to the 'Print' checkbox, 7 to the 'Printer' dropdown, 8 to the 'Email' checkbox, 9 to the 'Test Email' button, 10 to the 'From' input field, 11 to the 'Host' input field, 12 to the 'Password' input field, 13 to the 'Port' input field, 14 to the 'Colleague' section, 15 to the 'Laboratory' section, and 16 to the 'Host' input field.

No.	Description
1	Ticket Check this box if you want to print the ticket using the device's printer
2	header title If you'd like to include standard text at the top of each ticket, check the Header Title box, and then type your text in the lines below it. For example, you can enter your name and address on the lines.
3	Title text Enter title text here
4	Format Choose the type of printout you want
5	Symbol size Choose large or small print
6	print Check this box to print using an external printer. You can choose to print a Screenshot (Summary tab) and/or a Report.
7	Printer Select a default printer (only if you are connected to an external printer)
8	Email address Check the box to add email options.
9	Check email Select this button to send a test email. A confirmation message will appear
10	From Enter your email address in this field
11	Master Include the hostname here, for sample hostnames and configurations, go to #16, SMTP.

12	Password Enter this password in this field.
13	Port This is a field that is sometimes required according to your server configurations.
14	Colleagues Enter a colleague's name and email address here. Up to 3 colleagues can be included as the default recipient. To add more colleagues, press buttons 2 and 3.
15	Laboratory Enter the laboratory name and email address here. Up to 3 labs can be included as the default recipient. To add more labs, press buttons 2 and 3.
16	Host Configurations This popup shows possible host configurations such as SMTP configurations, port settings.


Email Address

To set up the email function, enter the following email addresses (#8), password (#12) and host (#11) and port (13). These specifications may vary depending on the specific email. Some common email specifications are listed here (#16) but may differ from those listed here. We recommend using the "Check Email" button to confirm that your email address is set up correctly. Also, make sure your device is connected to the internet.

Note: Emails cannot be sent automatically to ensure patient privacy. All emails can only be sent via the "Export" popup from the patient results.

Automatic actions

On the Automatic Actions tab of the Export tab, specify which export options will run automatically when this action is complete, and which options will be displayed in the Export dialog box.


configure automatic action settings:

- In the "Automatic actions" section (No. 1), select the options that should be performed automatically after the measurement is completed (Data export (No. 2) or Print (No. 3)).
- In the "Selected items for popup export" section (No. 4), select the options that should be displayed in the "Export" dialog box (options are available for exporting data (No. 5) or for printing (No. 6)).
- Select "Pop-up window always open" (7) to always open the "Export" dialog if "Export" is selected. Clear this option to automatically perform the actions selected in the selected Items to export pop-ups when export is selected.

Remote access

Remote access allows users on the same network to view results remotely using the device's IP address in a browser. In the browser, you can see a list of patients in chronological order, and for each patient, their test results. All security standards on your intranet may apply to the machine.

To select data to view on the remote access site:


No.	Description
1	Access to remote active When checked, remote access data will be generated automatically after each measurement. To remove this option, uncheck the box.
2	IP address Displays the IP address of the remote access data site.
3	Remote Access Data Recovery Click this button to add previously taken measurements to the remote access site.
4	Save Allows the user to save selected settings

If the Access Remote Active checkbox is checked, the data will be automatically exported after each measurement. It may take a few seconds for the dimension to appear on the remote access site. If you have data on the machine that is not yet on the remote site, or imported data, or old data that was not exported after measurement, you can export it retroactively.

There are two methods for exporting retro image data:

1. Export only missing measurements
2. Export the entire database.

Please note that option #2 will take some time and may slow down the device. Therefore, we recommend that these steps be performed while the device is idle.

To export retro image data, click the "Recover remote access data" button and the following screen will appear:

No.	Description
1	Cancel Deselects
2	Not Once selected, only the missing data will be added to the remote access site. Any data already stored on the remote access site will remain.
3	Yes This option will save the entire database on the remote access site, the new data as well as the previously existing data.


Note: Remote access (from intranet) only allows you to view images; changing data, such as replacing data, deleting data and making other changes is not possible remote access.

View in your browser

Once you've completed the setup, you'll see maps and data in your browser, or PC, MAC or tablet with any OS.

To view remote access in your browser:

- Open a browser and enter the IP address listed in the remote access menu.
- The following screen will appear:


No.	Description
1	IP address Enter the IP address listed in the remote access menu, may vary from time to time
2	A patient Patient name including birthday, gender and email address
3	Diagnostics A list of diagnostic measurements and dates was taken
4	IP address The identical address shown in the browser is an indication that you are properly connected to the site

5	Show name image Check the box to see the filenames of the images listed in the results.
6	Loading Select to download patient data

Please note that the database may take longer to load the first time.

- If you select the "show image names" option, you will see a red color on each map and data source file.
- Select the patient's name, select the measurement you want to load and press LOAD

In cases where images are not selected in the settings, they will not be saved and will be displayed NO DATA

- Press "EXIT" to return to the database.

Reports


The machine comes with a "Report" function that allows you to view the patient's test results clearly and concisely. You can choose the layout and data to be included in the report from over six predefined templates, or you can define your own custom template.

Report templates included with the machine:

- cataract
- aberrations
- topography
- glaucoma
- Mixed / all diagnostics

For more information about creating a report, see the "Creating a report" section.

The following is an example of a default report template:

No.	Description
1	Patient Information Includes patient name or ID, date and time of examination
2	Map Displays a map of the selected preference/default setting
3	Map name Displayed map name
4	Refraction Displays refraction data table
5	Sim K Displays Sim K Datasheet
6	Pachy Displays pachy datasheet: pachy thickness, idol angle values and front chamber depth
7	Tono Displays tono data table along with IOPc correction


Note: For a report to be available in the Export pop-up window on the Results screen, it must be selected in Settings > Export > Print: Report.

To learn how to define a custom report, go to Customizing a Report

Report setup

The machine offers the option to define a custom report. Three screens need to be configured to complete report customization: header and logo configuration, map configuration, and datasheet configuration.

To set up a report:
- In the "Settings" menu, select "Report":


No.	Description
1	Name Displays the name of the template shown in the window (#1)
2	Logo Choose to add a logo (jpeg or bmp file) from your network or USB
3	Template Name Add template name
4	impersonate Add text here for personalization
5	Add Select "Add" to add a new template, after you select this button, it will change to Cancel, and you can select it to cancel the template
6	Delete Choose to delete a template
7	Previous Select to display the previous screen
8	Next Select to proceed to the next step

- After using the "Add" button to open a new template, add a name, logo (optional) and text in fields 2, 3 and 4


- Select next

This screen allows you to customize the layout and card template:

No.	Description
1	Name Displays the name of the new template
2	Map information Select to change: type, eye selection and scale
3	Map Preview Shows a preview of the map, select the map you want to add to
4	Page layout Choose Landscape or Portrait
5	Images per page Choose the number of images per page, from 1 to 6 images
6	Page Size Choose your preferred page
7	Font size Choose your preferred font
8	Add Choose to add a new template
9	Delete Choose to delete a template
10	Previous Select to go to the previous screen
11	Next Select to proceed to the next step

- Select page layout settings in steps #6, 7 and 8
- To make changes to each map, you must first select it in window #3, then use the drop-down menus in #2 to change the map: size, map, type, and scale. Only the selected card will be changed.
- Select next
- The report data table screen allows you to make changes to the data tables:


No.	Description
1	Name Displays the name of the new template
2	Preview window Displays a preview of the report template.
3	Side Choose which eye to add the data table to, right or left
4	Map Select to select diagnostics to display in the datasheet: Topo, WF, Pachy and Tono
5	Type of Choose which type of datasheet will be displayed for each diagnostic, for example, if you select Topo (in #5), you can choose: Kerataconus or RMS data
6	Font size Choose the font size according to your preference.
7	Add, edit and remove Select the table you wish to edit and select 1 of these options according to your preference.
8	Add Template Select to add a report template.
9	Delete Template Select to delete the report template.
10	Previous Select to go to the previous screen
11	Save Select to save template
12	Cancel Select to undo changes

To add a data table to the template, select Add.

- In fields 3, 4 and 5 select the data that you would like to show in the table: displayed eye, diagnostics and data table type
- To make changes to existing tables, select the table in window #2 and then use the Edit/Delete buttons
- When you're d1 with your changes, select Save.

Note. Before creating a report, make sure it is selected in the Settings menu.

Once you have defined your own template according to your preferences, it will automatically appear in the report view menu shown here to the right of the dialog box below:

- Select the template of your choice from the menu, a preview of the report will appear in the window.
- Select "Export" to export to an external source or "Print" according to your preference.
- For more information about creating a report, go to the "Creating a Report" section

Restoring factory settings


You can return all configuration settings to their factory defaults.

To restore factory settings:

- On the "Configuration" screen, select "Factory". You will be asked to confirm that you want to restore all factory settings.
- Select Yes. Factory settings are restored for all configuration parameters.

Main10ance screen

The main10ance screen has many tabs that you can access to change settings in general device main10ance:

No.	Description
1	Help Choose to open H2O and get help
2	Version Selecting this button displays the software version number
3	Software update Selecting this option allows you to update the software, for more information go to Software Update
4	Save measurement files Saves last measurement for spy files
5	Topo Calibration testing Select to check topo calibration. This test must be carried out by an authorized technician.
6	Testing T1 calibration Artificial eye Select to check the tonal calibration of the artificial eye. This test should only be carried out by an authorized technician.
7	Topo Pachy Calibration Testing Select to check calibration for topo and pachy measurements. This test should only be conducted by an authorized professional.
8	Topo Pachy Calibration Testing Select to check calibration for topo and pachy measurements. This test should only be conducted by an authorized professional.
9	Advanced settings Select to open advanced settings, password required.
10	Engine reset Select to reset motors
11	Patient Directory Select to link to the patient in the database according to the first Title or ID
12	Database Path This shows the path to the database
13	Index setup This button updates the patient screen index. This is recommended after deleting/importing patient records
4	Data import

	Select this button to transfer data from USB or network
15	Search Choose to retrieve data from USB or network
16	Backup Select this button to backup the database
17	Deleting data Choose to delete your database.

Problems and solutions

You may encounter some of the following common problems while using the machine. If so, try the suggested solutions listed below. If the issue persists and the suggested solution does not resolve the issue, contact a qualified technical support representative or local distributor.

Problem	Possible reasons	Suggested Solutions
images are not correctly centered	- External lighting or windows near ism - Large spots light next to the block	- Make sure there are no bright lights or windows near the device or the patient. - Make sure there are no large irregular spots of light on the patient or on the device.
Measured pupil size different from expected or actual pupil size	- External lighting or windows near ism - Large spots light next to the block	- Make sure there are no bright lights or windows near the device or the patient. - Make sure there are no large irregular spots of light on the patient or on the device.
high order aberrations, which different from what expected	direct light incidence on the device or patient	- Make sure there are no bright lights or windows near the device or the patient. - Make sure there are no large irregular spots of light on the patient or on the device.
Tickets are not printed	- Paper jam - Out of paper	- Open the printing paper door and release paper. - Replace paper roll.

Service

This chapter explains how to clean the machine unit and how to perform current maintenance tasks.

-
- Device cleaning
 - Replacing printer paper

Cleaning the device

Important!

Before cleaning the device, turn it off and unplug it from the electrical outlet.

To clean the plastic surface of the machine, moisten a cloth with a commercial, non-abrasive cleaner and gently wipe the top, bottom, and front surfaces.

CAUTION: Do not spray or pour liquid directly onto the device.

CAUTION: Do not use caustic or abrasive cleaners

Chinrest

Regular cleaning of the chinrest with a soft cloth and alcohol is highly recommended. Use chin rest paper after each patient.

Headrest

Regular forehead cleaning with a soft cloth and alcohol is highly recommended. The forehead support is the only part of the machine that comes into contact with the patient.

Replacing Printer Paper






The printer is installed on the right side of the machine. If the paper roll, the printer's LED indicator flashes.


To insert a paper roll into the printer:

- Raise the handle in the middle of the paper compartment and pull the cover.
- If the printer has an empty paper roll, remove it.
- Insert a new roll with the paper end at the top of the roll.
- Load paper into the slot at the top of the compartment.
- Press the paper compartment door.

Main Information

- Screen contamination will affect measurement accuracy. Do not touch it with your hands and keep it away from dust.
- Remove smudges or dust from the screen with a soft cloth.
- Observe the operating and storage conditions.

Work		
Storage		
Transportation		

Do not install the device near TV or radio receivers as they may cause interference.

- If water or foreign objects get inside the appliance, unplug the power cord and contact your distributor. In the event of a malfunction (noise, smoke, etc.), turn off the power immediately and contact your distributor. Otherwise, fire or injury may result.
- Do not attempt to disassemble the instrument yourself, as this may cause a malfunction or fire. If operation fails, do not touch the inside of the appliance. Disconnect the power cord and contact your distributor.

EMC Compliance

The Comprehensive Eye Diagnosis Machine Models WAM700+ and WAC700 fully complies with relevant international and national laws.

When properly used, the Comprehensive Eye Diagnosis Machine models WAM700+ and WAC700 comply with international and national electromagnetic compatibility (EMC) laws and standards. These laws and standards set out the dosage of electromagnetic radiation allowed for the operation of the equipment and the ability of the equipment to withstand electromagnetic interference from the outside.

Equipment classification according to IEC 60601-1-2

Table 1 - Guidance and manufacturer's declaration - electromagnetic emissions

Guidance and manufacturer's declaration - electromagnetic emissions		
The Comprehensive Eye Diagnostic Machine Models WAM700+ and WAC700 is intended for use in the electromagnetic environment specified below. The purchaser or user of the Comprehensive Eye Tester Models WAM700+ and WAC700 should ensure that it is used in the specified electromagnetic environment.		
EMI test	Conformity	Electromagnetic environment - guidance
The group to which ME belongs according to CISPR 11 (GOST R 51318.11-2006)	Group 1	The Comprehensive Eye Diagnostic Machine Models WAM700+ and WAC700 uses RF energy for internal functions only. The level of radio frequency interference emissions is low and is not likely to cause malfunction of nearby electronic equipment. The Comprehensive Eye Diagnostic Machine Models WAM700+ and WAC700 is suitable for use in domestic premises and in premises directly connected to low-voltage electrical distribution networks (general-purpose electrical networks).
The class to which ME belongs according to CISPR 11 (GOST R 51318.11-2006)	Class B	
Harmonic current limits according to IEC 61000-3-2 (GOST 30804.3.2-2013)	Class A	
Voltage fluctuations and flicker according to IEC 61000-3-3 (GOST 30804.3.3-2013)	Corresponds	

Table 2. Guidance and manufacturer's declaration - immunity

Guidance and manufacturer's declaration - immunity			
The Comprehensive Eye Diagnostic Machine Models WAM700+ and WAC700 is intended for use in the electromagnetic environment specified below. To the purchaser or user The device for complex diagnostics of the eye, models WAM700+ and WAC700 ensure that it is used in the specified electromagnetic environment.			
Immunity test	Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharges (ESD) according to IEC 61000-4-2 (GOST 30804.4.2-2013)	±6 kV - contact discharge	±6 kV - contact discharge	The floor in the room is made of wood, concrete or ceramic tiles. With floors covered with synthetic material, the relative humidity of the air is at least 30%
	±8 kV - air discharge	±8 kV - air discharge	
Nanosecond impulse noise according to IEC 61000-4-4 (GOST 30804.4.4-2013)	±2 kV - for power lines	±2 kV - for power lines	Typical industrial or medical network.
	±1 kV - for input/output lines	-	
High energy microsecond impulse noise according to IEC 61000-4-5 (GOST R 51317.4.5-99)	±1 kV when interference is applied according to the "wire-wire" scheme	±1 kV when interference is applied according to the "wire-wire" scheme	For stable operation, it is necessary to use external surge protection devices, voltage stabilizers or uninterruptible power supplies!

	± 2 kV when interference is applied according to the "wire-to-ground" scheme	-	
Voltage dips, short interruptions and voltage changes in the input power lines according to IEC 61000-4-11 (GOST 30804.4.11-2013)	<5% U _n (voltage dip >95% U _n) for 0.5 period	<5% U _n (voltage dip >95% U _n) for 0.5 period	Typical industrial or medical network. For continuous operation under power surge conditions, it is recommended to use an uninterruptible power supply or a battery.
	40% U _n (voltage dip 60% U _n) for 5 periods	40% U _n (voltage dip 60% U _n) for 5 periods	
	70% U _n (voltage dip 30% U _n) for 25 periods	70% U _n (voltage dip 30% U _n) for 25 periods	
	<5% U _n (voltage dip >95% U _n) for 5 s	<5% U _n (voltage dip >95% U _n) for 5 s	
Lower frequency magnetic field (50/60 Hz) according to IEC 61000-4-8 (GOST R 50648-94)	3A/ m	3A/ m	Magnetic fields from the mains should not exceed typical levels for typical industrial networks.
NOTE: U _n - the voltage level of the electric network before the test action is applied.			

Table 3. Guidance and manufacturer's declaration - immunity

Guidance and manufacturer's declaration - Immunity NOT SUPPLIED FOR VITAL FUNCTIONS			
The Comprehensive Eye Diagnostic Machine Models WAM700+ and WAC700 is intended for use in the electromagnetic environment specified below. To the purchaser or user The device for complex diagnostics of the eye, models WAM700+ and WAC700 ensure that it is used in the specified electromagnetic environment.			
Equipment stability test	Challenge Level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF equipment, incl. cables, should not be used near the device closer than the recommended distance, calculated according to the formula according to the frequency of the transmitter:
Recommended Distance			
Conducted interference induced by radio frequency electromagnetic fields according to IEC 61000-4-6 (GOST R 51317.4.6-99)	3B from 150kHz to 80MHz	3B from 150kHz to 80MHz	$d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$

Radio frequency electromagnetic field according to IEC 61000-4-3 (GOST 30804.4.3-2013)	3V/m 80MHz to 2.5GHz	3V/m 80MHz to 2.5GHz	$d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz
			$d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.5 GHz

Where P is the maximum output power of the transmitter in watts (W) according to the manufacturer's specifications, and the recommended distance in meters (m).

d - recommended spacing, m;

Field strengths for propagation of radio waves from fixed radio transmitters, as observed in the electromagnetic environment ^{a)}, should be lower than the compliance level in each frequency band ^{b)}

Interference effects may occur in the vicinity of equipment marked with the symbol :

a) Field strengths for propagation of radio waves from fixed radio transmitters, such as base stations for radio teleph1 networks (cellular/wireless), and land mobile radios, amateur radio, AM and FM broadcast transmitters, and television transmitters cannot be determined by calculation with sufficient accuracy. To do this, practical measurements of the field strength must be made. If the measured values at the location of the Comprehensive Eye Diagnosis Unit , models WAM700+ and WAC700 exceed applicable compliance levels, the Comprehensive Eye Diagnostic Device , models WAM700+ and WAC700 should be observed to verify normal operation. If a deviation from normal functioning is detected during the observation, additional measures may need to be taken, such as reorientation or relocation of the Comprehensive Eye Apparatus , models WAM700+ and WAC700 .

b) Outside the frequency range from 150 kHz to 80 MHz, the field strength must not exceed (V 1) V/m.

Table 4. Recommended separation distances between portable and mobile RF communications equipment and the system

Recommended separation distances between portable and mobile RF communications equipment and the system			
NOT SUPPLIED FOR VITAL FUNCTIONS			
Device for complex eye diagnostics, models WAM700+ and WAC700 is included for use in an electromagnetic environment in which radiated interference levels are monitored. The customer or user of the system can help prevent electromagnetic interference by maintaining a minimum separation between portable and mobile RF communications equipment (transmitters) and the system as recommended below, taking into account the maximum output power of the communications equipment .			
Maximum transmitter output power (W)	Distance depending on transmitter frequency (m)		
	150 kHz ÷ 80 MHz d = $\left[\frac{3,5}{V_1} \right] \sqrt{P}$	80 MHz ÷ 800 MHz d = $\left[\frac{3,5}{E_1} \right] \sqrt{P}$	800 MHz ÷ 2.5 GHz d = $\left[\frac{7}{E_1} \right] \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74
1	1,17	1,17	2,33
10	3,69	3,69	7,38
100	11,67	11,67	23,33
When determining the recommended values of the spatial separation d for transmitters with a rated maximum output power not listed in the table, the expressions given are substituted with the rated maximum output power P in watts as specified in the transmitter manufacturer's documentation.			
NOTE <i>1 At 80 and 800 MHz, the higher field strength applies. 2 The expressions given may not apply in all cases. The propagation of electromagnetic waves is affected by absorption or reflection from structures, objects and people. 3 When determining the recommended separation distances d for transmitters with a rated maximum output power not listed in the table, the rated maximum output power P in watts specified in transmitter manufacturer's documentation</i>			

The Comprehensive Eye Diagnosis Machine Models WAM700+ and WAC700 requires special precautions regarding electromagnetic compatibility and must be installed and put into service in accordance with the EMC information contained in this document; Otherwise, the equipment may be adversely affected by mobile RF transceivers.

The use of accessories, transducers and cables other than those specified in this manual, with the exception of transducers and cables sold by the manufacturer as replacement parts for internal components, may increase radiation or reduce the life of the device.

The WAM700+ and WAC700 Comprehensive Eye Diagnosis Machine should not be used in close proximity to or in conjunction with other equipment.

Fuse replacement


Warning

To avoid electric shock, unplug the power cord from the unit before replacing the fuse.

Remove the holder with the blown fuse from the device by pushing it in and turning it counterclockwise .

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CONTENTS OF DELIVERY:

Device for complex eye diagnostics WAM700+:

- WAM700+ main unit – 1 pc.
- Network cable 220 V - 1 pc.
- Instruction manual - 1 pc.
- Dustproof cover - 1 pc.
- Chin rest paper - 1 set (1000 pcs in 1 set)
- Pin for the front stop - 2 pcs.
- Paper for the built-in thermal printer - 3 rolls.

Device for complex eye diagnostics WAC700:

- Main unit WAC 700 – 1 pc.
- Network cable 220 V - 1 pc.
- Instruction manual - 1 pc.
- Dustproof cover - 1 pc.
- Chin rest paper - 1 set (1000 pcs in 1 set)
- Pin for the front stop - 2 pcs.
- Paper for the built-in thermal printer - 3 rolls.

Packaging on the territory of the Russian Federation, in accordance with the standards of the Russian Federation:

The product must be packed and fixed in it during transportation.

Each package with the product, together with operational documentation, must be placed in a box in accordance with GOST 12301 or wrapped in paper in accordance with GOST 23436. The box or package must be tied with twine in accordance with GOST 17308 or pasted over with paper tape in accordance with GOST 18510, GOST 23436 or GOST 2228, tape paper-based adhesive according to GOST 18251 - or polyethylene tape with an adhesive layer according to GOST 20477 so that the package could not be opened without breaking the integrity.

For transportation, the preserved product in packaging with operational documentation must be placed in boxes made of sheet wood materials in accordance with GOST 5959 or plank. Plank boxes must be lined inside with waterproof material in accordance with GOST 2697.

Note:

1. Transportation of the product in railway containers can be carried out in primary packaging.
2. When shipping two or more transport places to 1 address, they can be packaged at the request of transport organizations. Maximum dimensions and w8 of the package according to GOST 24597.

Storage and transportation:

Transportation of the product can be carried out by all types of covered vehicles in accordance with the requirements of GOST R 50444 and the rules for the carriage of goods in force for each type of transport.

Transportation of cargo by sea should be carried out in accordance with the "Safety Rules for the Maritime Transportation of General Cargoes". Type of shipment - containers and small shipment.

The product, during transportation, must be resistant to climatic influences at an ambient temperature of -40 to + 70 C and humidity 10-95%.

The product, during storage, must be resistant to climatic influences at an ambient temperature of -10 to + 55 C and humidity 10-95%.

The product, during transportation, must be resistant to mechanical stress, with the following parameters:

- vibration loads - frequency range 10-55 Hz with a displacement amplitude of 0.35 mm;
- impact loads - peak impact acceleration 100 ms^{-2} with duration of impact acceleration 16 ms, impact frequency up to 120 min^{-1} .

Environmental requirements

The main type of possible hazardous impact on the environment is pollution of the atmospheric air of populated areas, soils and waters as a result of:

- emergency leaks (placers) of production materials;
- unorganized disposal of waste on the territory of the manufacturer or outside it;
- arbitrary dumping of waste in places not intended for this purpose.

The materials used in the manufacture of the product must not pose a danger to life, health of people and the environment both during operation and after the end of its service life.

manufacturer's warranty

The product is guaranteed for a period of 24 months from the date of sale from the manufacturer's warehouse. It is guaranteed that the product is free from defects related to production and the quality of production materials. Any defective part will be replaced if the product has not been tampered with and the product has been properly used during this period. Make sure that any malfunction is not due to

inadequate care of the product or failure to follow the instructions. In the event of a malfunction, subject to the operating requirements in accordance with the operating instructions, the manufacturer performs maintenance or replacement of the product. Maintenance of the product must be carried out by a specialist.

Official representative of the manufacturer: Essilor-LUYS-Optika Limited Liability Company, legal address: 127273, Moscow, st. Otradnaya, 2B, building 1, room No. 3.11, office. No. 51, fl. 3, tel/fax +7 (495) 137-48-48

Information about recycling.

Disposal of the product must be carried out in accordance with SanPiN 2.1.3684-21.

According to SanPiN 2.1.3684-21, the product belongs to class A - epidemiological safe waste.

The product must be disposed of if:

- end of service life;
- confirmation of facts and circumstances that endanger the life and health of health workers and indicate that the prescribed appointment has not been fulfilled.

Recycling is also subject to all packaging, including transport. Paper, wood, polyethylene and plastic must be disposed of separately in accordance with the rules for the collection, accounting and disposal of packaging.

Marking

The product must show:

- product name;
- designation of type and model;
- manufacturer;
- manufacturer's address;
- supply voltage;
- frequency of the supply voltage;
- power consumption;
- electrical safety classification;
- serial number;
- Symbols

Symbol	Description
	Important: refer to the instruction manual
	Alternating current
	Type BF
	The equipment must be returned to the manufacturer or disposed of
	Manufacturer
	Production year
	Serial number
	Compliance with medical directive 93/42/EC

The packaging must indicate:

- series;
- year of issue (last two digits)
- trademark of the manufacturer;
- product name;
- year and month of packing;
- warning sign: "Caution! Fragile!"

List of applicable safety and electromagnetic compatibility standards

GOST R 50444-92 Medical instruments, apparatus and equipment. General specifications

GOST IEC GOST R IEC 60601-1-2010 Medical electrical products. Part 1. General safety requirements, taking into account the main functional characteristics.

GOST 31590.1-2012 Ophthalmic instruments. Part 1: General requirements for ophthalmic instruments and test methods

GOST ISO 14971-2011 Medical products. Application of risk management to medical devices

GOST R IEC 62304-2013 Medical products. Software. Life cycle processes

GOST R IEC 60601-1-2-2014 Medical electrical products. Part 1-2. General safety requirements, taking into account the main functional characteristics. Parallel standard

GOST 28312-89 Professional radio-electronic equipment. Conditional graphics

GOST R IEC 60601-1-8-2007 Medical electrical products. Part 1-8. General safety requirements. General requirements, tests and guidelines for the use of alarm systems for medical electrical devices and medical electrical systems

GOST 30804.3.2-2013 (IEC 61000-3-2:2009) Electromagnetic compatibility of technical means. Emission of harmonic components of the current by technical means with a current consumption of not more than 16 A (in 1 phase). Norms and test methods

GOST 30804.3.3-2013 (IEC 61000-3-3:2008) Electromagnetic compatibility of technical means. Limiting voltage variations, voltage fluctuations and flicker in general low voltage power supply systems. Technical means with a current consumption of not more than 16 A (in 1 phase), connected to the electrical network if certain connection conditions are not met. Norms and test methods

GOST 30804.4.2-2013 (IEC 61000-4-2:2008) Electromagnetic compatibility of technical means. Resistance to electrostatic discharges. Requirements and test methods

GOST 30804.4.3-2013 (IEC 61000-4-3:2006) Electromagnetic compatibility of technical means. Resistance to radio frequency electromagnetic field. Requirements and test methods

- GOST 30804.4.4-2013 (IEC 61000-4-4:2004) Electromagnetic compatibility of technical means. Resistance to nanosecond impulse noise. Requirements and test methods
- GOST R 51317.4.5-99 (IEC 61000-4-5-95) Electromagnetic compatibility of technical means. Immunity to microsecond impulse noise of high energy. Requirements and test methods
- GOST P 51317.4.6-99 (IEC 61000-4-6-96) Electromagnetic compatibility of technical means. Immunity to conducted interference induced by radio frequency electromagnetic fields. Requirements and test methods
- GOST R 50648-94 (IEC 1000-4-8-93) Electromagnetic compatibility of technical means. Resistance to the magnetic field of industrial frequency. Technical requirements and test methods
- GOST 30804.4.11-2013 (IEC 61000-4-11:2004) Electromagnetic compatibility of technical means. Resistant to dips, short interruptions and power supply voltage fluctuations. Requirements and test methods
- GOST R 51318.11-2006 (CISPR 11:2004) Electromagnetic compatibility of technical means. Industrial, scientific, medical and household (PNMB) high-frequency devices. Industrial radio interference. Norms and methods of measurements
- GOST 30805.14.1-2013 (CISPR 14-1:2005) Electromagnetic compatibility of technical means. Household appliances, electrical tools and similar devices. Industrial radio interference. Norms and methods of measurements
- GOST R 51318.15-99 (CISPR 15-96) Electromagnetic compatibility of technical means. Industrial radio interference from electrical lighting and similar equipment. Norms and test methods
- GOST 30805.16.1.2-2013 (CISPR 16-1-2:2006) Electromagnetic compatibility of technical means. Requirements for equipment for measuring the parameters of industrial radio interference and noise immunity and measurement methods. Part 1-2. Equipment for measuring the parameters of industrial radio interference and noise immunity. Devices for measuring radio interference and testing for immunity to radio interference
- GOST 30805.22-2013 (CISPR 22:2006) Electromagnetic compatibility of technical means. Information technology equipment. Industrial radio interference. Norms and methods of measurements

WARRANTY CARD

Product category _____
 Manufacturer _____
 Model _____ Stamp _____
 Serial number/art. _____ about sale _____
 Date of sale _____
 Organization of the seller _____
 Seller's signature _____
 The product has been tested and has no damage.
 I have read and agree to the terms of use.
 Buyer's signature _____
 The warranty period is 24 months from the date of sale from the manufacturer's warehouse.
 For warranty questions please contact:

Essilor-LUYS-Optika Limited Liability Company, legal entity address: 127273, Moscow, st. Otradnaya, 2B, building 1, room No. 3.11, office. No. 51, fl. 3, tel/fax +7 (495) 137-48-48

ATTENTION! Make sure that the warranty card is filled out correctly. The serial number of the product / article must exactly match the number in the coupon. All columns of the coupon (with the exception of the columns for the guarantee center) must be filled in.

In case of inaccurate or incomplete filling of the coupon, the warranty for the product is considered lost. The manufacturer's warranty is carried out in accordance with the Civil Code of the Russian Federation on a general basis.

The product must not bear traces of mechanical or thermal effects, as well as traces of exposure to various chemical agents. The product may only be operated in accordance with the instruction manual for the product.

WARRANTY STATEMENT

THE DATE	DESCRIPTION	WARRANTY CENTER

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