
MILESTONE
H E L P I N G
P A T I E N T S

Approved

(Diego Cortesi - CEO)

MILESTONE S.r.l.
Via Falebenefratelli, 1/5
24010 SORISOLE (BG) - Italia
P. IVA n. C. F. 01679330163

SealSAFE

Automatic Fixative Filling and
Vacuum Sealing System in
models, with accessories

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We certify that the signer is a duly authorized signatory on behalf of the Company for any documents and acts used in trade with foreign countries.

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Monica Loisi

BERGAMO, 30/11/2021

Thank you for having selected our system and welcome to the ever growing world club of users for Milestone laboratory instrumentation.

We are sure that you will be completely satisfied with this new tool entering your laboratory.

We invite you to read carefully this operator manual and to keep it in reach for convenient and fast consulting.

For any possible clarification or any request for assistance please contact either our Representative in your country:

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Corso Europa Unita, 6

IT ~ 24030 Valbrembo (BG), Italy

Please read the user manual before using the device.

MM095 – SealSAFE – Operator Manual

SUMMARY

1. INTRODUCTION	5
1.1. Symbols used in this manual	5
1.2. Description and Designated uses	6
1.2.1. Designated uses	6
1.2.2. Delivery set	7
1.3. Technical specifications	8
1.3.1. SealSAFE structure	9
1.3.2. Classification	11
Classification as per GOST 20790-93:	11
Class depending on the possible consequences of failure during use is a class B device.	11
1.4. Touch control terminal	12
1.5. Transportation and storage conditions	12
1.6. Environmental conditions	13
1.7. Warnings and Precautions	13
1.8. Instrument marking	16
1.8.1. Package marking (wooden box)	16
1.8.2. Nameplate	18
1.8.3. Barcode reader marking	20
1.8.4. Location of symbols on the SealSAFE	22
1.9. Compatible reagents for the fixation	22
1.10. Barcode reader	23
2. INSTALLATION AND SETTING UP THE INSTRUMENT	24
2.1. Site requirements	24
2.2. Unpackaging and check list	25
2.2.1. Accessories	26
2.3. How to install the bag holder	26
2.4. How to install the barcode reader	27
2.5. Fume extraction system	29
2.5.1. How to install the External exhaust FORMALIN fumes filter kit	30
2.5.2. How to install the External Self Standing XL exhaust FORMALIN fumes filter kit	32
2.6. Electrical installation of SealSAFE	34
2.6.1. Power supply	34
2.6.2. Back side connection plate	34
2.7. Final placement	34
2.8. Safety Switch ON and OFF	36
2.9. Reagents	37
3. OPERATING SEALSAFE	38
3.1. Log IN and Log OUT	38
3.1.1. Create a new user	38
3.2. Set Date and Time	40
3.3. How to load SealSAFE with reagents	41
3.3.1. Reagents set up	41
3.3.2. How to replace an empty tank	44
3.4. Owner Settings	44
3.4.1. General	45
3.4.2. Ventilations	45
3.4.3. HEPA filter	45
3.4.4. FORMALIN filter	46
3.4.6. Enables	49
3.4.7. Milewatch	50
3.4.8. Print	50
3.4.9. Case ID	50
3.5. Run a program	52
3.6. Vacuum Only program	53
3.7. Vacuum + Fixative program	64
3.7.1. Reagent management during a "Vacuum + Fixative" program	76
3.7.2. How to use a vacuum bag more times for the same specimen	77
3.8. Power failure during a program	80
3.9. Retrieving saved processing results	80
3.9.1. How to delete a saved processing	82
3.9.2. Export processing data from the terminal	82

3.9.3. Install and operate with the LogViewer.....	84
3.10. How to modify a program.....	87
3.11. How to set up a Favorite program.....	89
3.12. Temperature Monitoring System.....	89
3.12.1. Sender Mode for SealSAFE - Sending specimens.....	90
3.12.2. How to enter Case IDs manually.....	91
3.12.3. How to enter Case IDs automatically.....	92
3.12.4. How to activate a card.....	94
3.12.5. How to activate a card in a Delay mode.....	95
3.12.6. Receiver Mode - Receiving specimens.....	96
4. ACCESSORIES.....	100
4.1. External exhaust FORMALIN fumes filter kit (CODE 61550/F).....	100
4.1.1. Handling and storage.....	100
4.1.2. How to change the filter.....	101
4.2. External Self Standing XL exhaust FORMALIN fumes filter kit (CODE 68620).....	102
4.2.1. Handling and storage.....	102
4.2.2. How to change the filter.....	103
4.3 HEPA filter (CODE 68290).....	104
5. CONSUMABLES.....	106
5.1. Vacuum bags.....	106
5.1.1. Storage of Vacuum bags.....	106
5.1.2. Stand up vacuum bags.....	106
5.1.3. Vacuum bags not Stand up.....	107
6. ADDITIONAL DEVICES THAT ARE NOT PART OF THE SYSTEM, BUT WHICH CAN BE USED WITH IT.....	109
6.1. Kit Refrigerator box.....	109
6.2. Vacuum bags.....	109
6.3. Temperature Monitoring System.....	109
6.4. Compatible reagents for the fixation.....	109
7. MAINTENANCE.....	110
7.1. Operating room personnel.....	110
7.1.1. Ordinary maintenance.....	110
7.1.2. Extraordinary maintenance.....	111
7.2. Service personnel.....	113
7.3. Calibration of the instrument.....	113
8. Responsibility of the Manufacturer.....	115
9. Repair and service information.....	115
10. WARRANTY.....	116
11. WASTE DISPOSAL.....	116
11.1. Waste disposal of Vacuum bags.....	116
11.2. Waste disposal of tanks.....	116
11.3. Waste disposal of fixatives.....	116
11.4. Waste disposal of the equipment.....	117
11.5. Waste disposal of filters (if present).....	117
12. DECONTAMINATION AND CLEANING.....	118
12.1. Safety precautions.....	118
12.2. SealSAFE decontamination and cleaning procedure.....	119
12.2.1. Detergents.....	119
12.2.2. Procedure.....	119
13. PRECAUTIONS IN CASE OF THE MEDICAL DEVICE MALFUNCTION, FAILURE, OR DEVIATIONS IN FUNCTIONING.....	122
A. APPENDIX A – “VACUUM+FIXATIVE” PROGRAMS GUIDELINES.....	123
A.1. How to correctly position a bag into the cavity.....	123
A.2. How to correctly seal a vacuum bag.....	124
B. APPENDIX B - “VACUUM ONLY” PROGRAMS SELECTION GUIDELINES.....	125
C. APPENDIX C – SPECIMENS TRANSPORTATION GUIDELINES.....	126
D. APPENDIX D – PRINTER AND BARCODE READER.....	127
D.1. How to use the printer.....	127
D.1.1. Storage of paper rolls.....	128
D.1.2. Changing the paper roll.....	128
D.1.3. Printer troubleshooting.....	129
E. APPENDIX E – BAG SEALER WITH HANDLE.....	130

F. APPENDIX F - RAPID CONNECTOR WITH HOSE SUITABLE FOR POLITAINER TANKS (CODE 68630)	13
G. APPENDIX G – THE LIST OF NATIONAL STANDARDS APPLIED BY THE MANUFACTURER.	134

1. INTRODUCTION

1.1. Symbols used in this manual

	Please read carefully this instruction.
	An instruction accompanied by this symbol provides a cautionary statement: WARNING . Failure to follow the instruction may endanger the user or cause damage to the instrument.
	An instruction accompanied by this symbol provides a caution against electric shock. Be sure to follow the instruction to prevent accidents resulting from electric shock.
	Biohazard.
	Surfaces identified with this symbol are at danger temperature (more than 70°C).
	Warning: Laser light
	CE logo: This instrument complies with European community directives.
	WEEE European directive symbol (2002/96/EC): Electric/electronic device do not throw out in the environment.
	IVD medical device according to 98/79/EC directive.
	Indicates a medical device for one-time use or for use upon an individual patient during a single procedure.
	Indicates the manufacturer's batch number.
	Indicates the manufacturer's serial number.
	Indicates the manufacturer's catalogue code.
	Switched on
	Switched off
	Alternating current
	Nominal voltage, V
	Maximum nominal power, W
	USB slot
	Manufacturer
	Fragile

	Used by... Indicates the date after which the device should not be used.
	Norma Oficial Mexicana (NOM) Certificate
	Electrical Safety Certification (indication of Mexican agency)
	Canadian Standards Association (CSA)
	KC (certification in Korea)
	China RoHS China RoHS environmental safety symbol. The number in the symbol means the "Environmental Friendly Usage Period" of the device in years.
	Indian standard
	Australia and New Zealand RCM (Regulatory Compliance Mark)
	Conformity mark to Eurasian standards
	Japan VCCI Certificate
	RoHS Mark (Chinese standard)
	UL Mark for North America

1.2. Description and Designated uses

The medical product name is "SealSAFE Automatic Fixative Filling and Vacuum Sealing System in models, with accessories" (hereinafter –"Automatic Fixative Filling and Vacuum Sealing System", "the SealSAFE", "the medical product", "the medical device", "system" and "instrument").

1.2.1. Designated uses

SealSAFE is an auxiliary equipment for in vitro diagnostics, specially designed and intended for packing the histological specimens of biological tissues and human organs into vacuum bags for subsequent examination in histopathological laboratories of medical institutions:

- Preservation and transportation of biological samples in a container or without it, fresh or fixed in a vacuum bag;
- Fixation before grossing with a defined quantity of fixative into the vacuum bag;
- Long term specimen preservation (if required with a defined quantity of fixative into the vacuum bag).

Description of a target analyte	N.A.
Functional purpose	Auxiliary equipment for in vitro diagnostics
Specific pathology, condition or risk factor for the detection, identification or differentiation of which the IVD device is intended	N.A.
Type of the specimen analyzed	Biological specimens of human organs and tissues

Population, demographic aspects of the medical device use	SealSAFE is not intended for a specific patient group. It is allowed to use the device for biological samples of human organs and tissues of all populations, patient groups of different age.
---	--

Contraindications: not applicable

Side effects: not applicable

Type of samples intended for sealing by the SealSAFE: any biological specimens of human organs and tissues

The SealSAFE may be operated by trained laboratory personnel only (clinical laboratory specialist or assistant). All laboratory personnel designated to operate this instrument must read the Operator Manual carefully and must be familiar with all technical features of the instrument before attempting to operate it.

Only for the professional users.

The instrument is intended for laboratory use (indoor use) only.

For clinical pathology tests.

The Milestone SealSAFE must be operated exclusively with the reagents listed in section "Compatible reagents for the fixation" 1.9.

The Milestone SealSAFE must be operated exclusively as reported in this manual.

Any other use of the instrument is considered improper and may void the warranty provided by the manufacturer.

1.2.2. Delivery set

SealSAFE Automatic Fixative Filling and Vacuum Sealing System in models, with accessories

Models:

- I. SealSAFE Automatic Fixative Filling and Vacuum Sealing System, basic;

Composed of:

- 1) SealSAFE Automatic Fixative Filling and Vacuum Sealing System, basic - 1 pc.;
- 2) Exhaust Tube complete - 1 pc.;
- 3) Metal clamp - 2 pcs;
- 4) Main power cable - 1 pc.;
- 5) Barcode reader mounted on its support - 1 pc.;
- 6) Bag holder - 1 pc.;
- 7) USB with "LogViewer" software included - 1 pc.;
- 8) Pipes numbered: 1, 2, EXT1 and EXT2 - 4 pcs;
- 9) Plastic rings 1-2 - 2 pcs;
- 10) 3 Different sizes of cups for the tanks of the drawer (small size diameter 28-33 mm, medium size diameter 32-38 mm, big size diameter 42-49 mm) - 6 pcs;
- 11) Metal rings - 2 pcs;
- 12) Allen wrench - 1 pc.;
- 13) Metallic plate - 1 pc.;
- 14) Rapid connector with hose suitable for Politainer tanks (if required) - 1 pc.;
- 15) Rolls of papers for the printer (if required) - 10 pcs;
- 16) Plastic stylus to work with the touch-screen terminal - 1 pc.;
- 17) Operator manual in electronic format - 1 pc.;
- 18) Bag sealer with handle (if required) - up to 10 pcs;
- 19) Operator manual in paper format - 1 pc.;

II. SealSAFE WDLR Automatic Fixative Filling and Vacuum Sealing System, with integrated slot for temperature monitoring system;

Composed of:

- 1) SealSAFE WDLR Automatic Fixative Filling and Vacuum Sealing System, with integrated slot for temperature monitoring system - 1 pc.;
- 2) Exhaust Tube complete - 1 pc.;
- 3) Metal clamp - 2 pcs;
- 4) Main power cable - 1 pc.;
- 5) Barcode reader mounted on its support - 1 pc.;
- 6) Bag holder - 1 pc.;
- 7) USB with "LogViewer" software included - 1 pc.;
- 8) Pipes numbered: 1, 2, EXT1 and EXT2 - 4 pcs;
- 9) Plastic rings 1-2 - 2 pcs;
- 10) 3 Different sizes of cups for the tanks of the drawer (small size diameter 28-33 mm, medium size diameter 32-38 mm, big size diameter 42-49 mm) - 6 pcs;
- 11) Metal rings - 2 pcs;
- 12) Allen wrench - 1 pc.;
- 13) Metallic plate - 1 pc.;
- 14) Rapid connector with hose suitable for Politainer tanks (if required) - 1 pc.;
- 15) Rolls of papers for the printer (if required) - 10 pcs;
- 16) Plastic stylus to work with the touch-screen terminal - 1 pc.;
- 17) Operator manual in electronic format - 1 pc.;
- 18) Bag sealer with handle (if required) - up to 10 pcs;
- 19) Operator manual in paper format - 1 pc.;

Accessories:

- 1) External exhaust FORMALIN fumes filter kit - 1 pc.;
composed of:
 - 1.1) Filter housing - 1 pc.,
 - 1.2) Filter - 1 pc.,
 - 1.3) Label to the filter - 1 pc.,
 - 1.4) Kit of 3 different lengths exhaust tubes complete with metal clamps - 1 pc.,
- 2) External Self Standing XL exhaust FORMALIN fumes filter kit - 1 pc.;
composed of:
 - 2.1) Filter housing - 1 pc.,
 - 2.2) exhaust FORMALIN filter XL - 2 pcs,
 - 2.3) Label to the filter - 1 pc.,
 - 2.4) Filter dividing chamber - 1 pc.,
 - 2.5) Exhaust tube complete with metal clamps - 1 pc.
- 3) HEPA filter - 1 pc.

1.3. Technical specifications

Model: SealSAFE Automatic Fixative Filling and Vacuum Sealing System, basic:


Code 68460: Power supply input: 230V±10% ~ 50Hz 1500W, 1800Wp

Model: SealSAFE WDLR Automatic Fixative Filling and Vacuum Sealing System, with integrated slot for temperature monitoring system (with Temperature data logger):

Code 68463: Power supply input: 230V±10% ~ 50Hz 1500W, 1800Wp


1.3.1. SealSAFE structure

General view:

INTERNAL CAVITY

Opening the main transparent cover you can identify the following components:

**Legend:**

- 1- Silicon Bar
- 2- EPDM seal
- 3- Hole formalin discharge
- 4- Vacuum pump switch
- 5- Sidebar and metal plate for positioning bags
- 6- Extract air grid
- 7- Load Cell
- 8- Welding bar
- 9- Nozzle
- 10- Grid for cavity ventilation
- 11- Terminal
- 12- Printer

Back side connection plate

- 1) Ethernet / LAN Plug
- 2) Main power supply
- 3) EXT1 outlet external fixative source 1
- 4) EXT2 outlet external fixative source 2

1.3.2. Classification

Potential risk class for SealSAFE - 2a

Type of the medical device in accordance with the Russian nomenclature classification is 116950.

Classification as per GOST 20790-93:

Class depending on the possible consequences of failure during use is a class B device.

Group depending on mechanical stress - device of group 2

IP protection degree (according to GOST 14254) - IPX0.

Laser product class (barcode reader) - 2.

Technical parameters:

Operating programs of instrument:	2 types of programs: - Vacuum Only (VO) - Vacuum+Fixative (VF)
Max. acceptable time for setting the operating mode:	up to 30 s
Average time of sample evacuation for the program: Vacuum Only (VO) P1 P2 P3 P4	Regardless of the size of the vacuum bag used and the biological sample being sealed P1 - 40 s ± 5 s P2 - 33 s ± 5 s P3 - 26 s ± 5 s P4 - 13 s ± 5 s
Average time of reagent filling and sample evacuation for the program: Vacuum+Fixative (VF)	1 min 25 s ± 5 s
Accuracy of weighing and fixative dispensing (for all ratios of specimen weight: fixative weight)	± 10 g
Force applied when blocking the wheels of the device	Not more than 35 N
Force applied when moving the SealSAFE	Not more than 40 N
Installed vacuum level:	Installed VO vacuum level: P1 – 25 mbar P2 – 50 mbar P3 – 75 mbar P4 – 325 mbar
Quantity of tanks:	• 2 tanks with a capacity up to 20L maximum and dimensions (LxHxW) not more than (73,5 x 50,5 x 42) cm can be placed in the drawer, • 2 external tanks/containers of the desired capacity and dimensions can be placed next to the instrument.
HEPA filter:	Filter elements of the "BSV filtration class DA" series are classified HEPA ISO 35H. Retention efficiency: >99.999 % (produce a sterile effluent)
External exhaust FORMALIN fumes filter:	Efficiency: 50% Formalin flow rate upstream of the filter 1.77 mg/m ³ – sampling test 3 L/h of formalin Formalin flow rate downstream of the filter 0.9 mg/m ³ sampling test 3 L/h of formalin Breakthrough: 300L of 10% formalin
External Self Standing XL exhaust FORMALIN fumes filter:	Efficiency: 90% Formalin flow rate upstream of the XL filter 1.77 mg/m ³ – sampling test 3 L/h of formalin Formalin flow rate downstream of the XL filter 0.17 mg/m ³ sampling test 3 L/h of formalin Breakthrough: 1000L of 10% formalin
The sound power level of noise:	Not more than 65 dB

Fan air flow range:	47-50 m³/h (cube meter per hour).
Fuses:	F1-F2: 8A 10x38 GG F7: 20A 10x38 AM
Barcode reader:	Emitted wavelength - 650 nm; output radiation - 1 mW.
Characteristics of the seam of the vacuum bag	The seam should be even, uniform, without gaps, burnt spots and folds.


The Quality System Management of Milestone s.r.l. is ISO 9001:2015 and ISO 13485:2016 certified.

1.4. Touch control terminal

The following picture shows the Touch Control TERMINAL of the SealSAFE unit. It is possible to operate the unit by simply touching the icons on the screen.

Two USB ports are located on the lateral sides of the Touch Control Terminal, one on the right and one on the left (see the arrow in picture below, right side).

When you don't use the USB ports, it's required to leave the USB protection inserted into the slot, as shown in the picture above.

Characteristics:

TFT-LCD 7" widescreen color display
Resolution 800X480,
Resistive Touchscreen
Mass storage 4GB
1CPU 1GHz

Software:

Windows™ Embedded Compact 7

Programs:

Multiple pre-stored programs, with customizable user defined programs.

I/O-Interface:

2x USB port 2.0, 1x LAN Ethernet 10/100 Mbps RJ-45.

Software safety class in accordance with the standard IEC 62304-2006: Class A: No injury or damage to health is possible.

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1.5. Transportation and storage conditions

The instrument can be transported by rail, road, river and air transport in covered vehicles.

Transportation and storage conditions:

Temperature: -20°C to +60°C

Humidity: up to 50% (at 40°C) non-condensing
Special requirements of the manufacturer.

Fastening of a transport packaging in vehicles and transportation of the products must be carried out in accordance with the rules applicable for the respective means of transport.

The arrangement and fastening of the transport boxes should ensure their stable position, and displacement during transportation is not allowed.

Since the instrument is large-sized and heavy, it is recommended to place it on the lower rows of racks in the manufacturer's (consumer's) warehouses.

The instrument must only be transported in an upright position.

Before turning SealSAFE on (when coming from a storage room), let the unit reach working environmental conditions (at least half an hour).

Rolls of papers for the printer (component)

Transportation conditions
Joint transportation with the instrument is allowed
Temperature: -20°C to +60°C
Humidity: up to 50% (at 40°C) non-condensing

Storage conditions
Store separately from the instrument
Temperature: +20 to 40°C.
Humidity: up to 50% (at 40°C) non-condensing.

1.6. Environmental conditions

SealSAFE should be operated indoors only.

Working temperature:	15°C - 30°C
Humidity (relative):	up to 80%
Maximum altitude:	2000m
Special requirements of the manufacturer.	

1.7. Warnings and Precautions

Should the equipment be used in different ways than those specified by the manufacturer, the warranty provided by the manufacturer may be void.

Main power plugs are considered disconnecting device. Disconnect all plugs from the socket outlet before assembling of the system, before connection of the accessories (if present) and before cleaning.

The instrument must not be placed in such a way that it is difficult to access the disconnecting device. Equipment shall be connected to a main power socket outlet which is grounded, accessible and visible by the operator. DO NOT use extension cords without a protective conductor.

Before starting work, while the device is in a safe state, it is necessary to visually make sure that the device is evenly placed on a horizontal surface, and check the absence of external damage to the device and power cable.

Equipment is provided with replaceable fuses, however their replacement is to be performed by properly authorized service personnel. If the system is not properly working, please contact the local authorized representative of the manufacturer.

This equipment has been designed and tested in accordance with the requirements of CISPR 11 (class A). In residential areas, it can cause radio interference, and in this case you should take measures to reduce the level of interference.

This equipment complies with the requirements for immunity and electromagnetic emission of GOST R IEC 61326-1-2014 and GOST R IEC 61326-2-6-2014, but when the equipment is connected to other objects or in other operating conditions, electromagnetic emission that exceeds the level specified in this standard may occur.

It is recommended to place the device in places where the intensity of external electromagnetic radiation and electrostatic discharge is below the maximum allowable intensity.

Do not use this equipment near sources of strong electromagnetic radiation, ionizing radiation, or electrostatic discharges, as they may disrupt its proper operation.

When connecting products or equipment that are not included in the delivery to the device, electromagnetic emission that exceeds the permissible level may occur.

It is recommended to use the device in the specified environmental conditions (see Chapter 1.6) in order to avoid adverse effects on its operation.

Filters and reagents after their use should be treated as a biohazardous material. They should be disposed of in accordance with federal, municipal and local government regulations (see chapters 11.3 and 11.5).

All parts of the equipment and all accessories must be supplied only by the manufacturer.

DO NOT OPEN THE CHASSIS. Failure to follow the instruction in the documentation might lead to a reduction of device safety. If the equipment is not working properly, please contact your supplier or the manufacturer.

Do not move the instrument, accessories, panels or cover. Only authorized and qualified **SERVICE** personnel may repair the instrument and access the instrument's internal components.

It is important that normal standards of safety and good laboratory practices are employed. Always use common sense and the best known practice when operating the instrument. The institution owning the instrument has primary responsibility for accident-free operation together with designated personnel who operate, service or repair it.

To avoid damage to the instrument, use only the reagents listed in section 1.9 Compatible reagents for the fixation. Failure to follow the instruction in the documentation might lead to a reduction of device safety.

Be careful: reagents exhaust are toxic and carcinogenic. The use of SealSAFE as described in this operator manual keeps the operator in safety. In case of accidental spillage of fixative, follow the local regulations that contain important information, as the use of suitable Personal Protective Equipment (PPE).

To avoid damage to the instrument, use only the vacuum bags described in this manual, supplied by Milestone s.r.l. For details, see the chapter 5.1. For disposal, see the chapter 11.1.

To prevent triggering of thermal protection in certain components, it is not recommended to run more than 3 cycles for 15 minutes for small samples.

During normal use, some surfaces of the equipment become hot. Contact with such surfaces without proper precautions can cause burns. At the end of a program, the sealing bar can be very hot.

Do not obstruct air discharge.

It is unacceptable to replace the detachable power cable with cables with inappropriate nominal characteristics.

During the operation of the device, its cover must be closed. Milestone S.r.l. is not responsible for any degradation in system performance caused by opening the cover during the process of packaging and sealing biological samples in vacuum bags. **CAUTION** when closing the cover: Risk of crushing! Do not put your hands in the cover movement area.


When using equipment with violations of the requirements established by the manufacturer, the level of equipment protection may be reduced.

Do not direct the barcode reader into your eyes. Laser light, do not stare into beam.
CLASS 2 LASER EQUIPMENT.

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

1.8.3. Barcode reader marking


For symbols explanation see Chapter 1.1.

1.8.4. Location of symbols on the SealSAFE


Back of the instrument:

Side of the instrument:

Cover of the instrument:

For symbols explanation see Chapter 1.1.

1.9. Compatible reagents for the fixation

1. Formalin 10%, buffered (Registration Certificate in the Russian Federation No. P3H 2013/1170, cl. 14 or other reagents with similar composition registered in the territory of Russia).
2. Formalin equivalent: "GreenFix" solution (Registration Certificate in the Russian Federation No. P3H 2013/1025, cl. 22 or other reagents with similar composition registered in the territory of Russia).

DO NOT USE REAGENTS DIFFERENT FROM THOSE LISTED ABOVE.

In case of doubt please contact:
application@milestonemedarl.com <mailto:application@milestonemedarl.com>

Reagents other than those listed above may damage some components of the Instrument. Check with authorized representative in your country the reagents allowed for circulation and their availability.


1.10. Barcode reader

Barcode reader is included in the delivery set of the unit.
It ensures an intuitive and fast reading of barcodes. Ergonomically designed for hands-free or handheld operations.

Refer to the chapter 2.4 to get the instructions on how to install it on the unit.

2. INSTALLATION AND SETTING UP THE INSTRUMENT

The SealSAFE is a precision instrument that requires the utmost care when it is unpacked and installed. Move SealSAFE to its final location. Make sure that the floor is level and is made of non-flammable material. Installation and setting up of the instrument should only be carried out by a specially trained engineer. The following pictures show the dimensions of SealSAFE and the space required.

Dimensions of the instrument with drawer and cover opened:

W x H x D: 55 cm x 108 cm (140 cm cover opened) x 85 cm (142 cm drawer opened) (± 2 cm)

After installing the device in its final location, it is not recommended to move it.

2.1. Site requirements

SealSAFE is delivered in a wood box with these dimensions:

Width: 77cm ± 5 %	Height: 134cm ± 5 %	Depth: 102cm ± 5 %
-----------------------	-------------------------	------------------------

The space dimensions required to locate the unit are the following (drawer and cover opened):

Width: 95cm ± 5 %	Height: 140cm ± 5 %	Depth: 185cm ± 5 %
-----------------------	-------------------------	------------------------

SealSAFE is a heavy instrument. Contact a structural engineer before placing the unit (analyze the building structure).

SealSAFE Weight (read carefully):

SealSAFE FULL (all tanks): 200kg ± 2 kg (2 tanks of 20L).


SealSAFE EMPTY: 160kg ± 2 kg (no tanks inside)

Unit density (when full): 425kg/m² ± 5 %.

2.2. Unpackaging and check list

When unpacking, check that all parts correspond to the packing list included in the shipment. Utmost care must be made during unpacking of SealSAFE, to avoid any scratch or damage to the external coating. The package includes the following items:

- Qty.1 Exhaust Tube complete (CODE 61377) and q.ty 2 metal clamp (CODE 50051),

- Qty.1 Main power cable (CODE 50036 for 230V version),
- Qty.1 barcode reader mounted on its support (to be fixed on the unit, see the chapter 2.4) (CODE 68589A)

- Qty.1 bag holder (to be fixed on the unit, see the chapter 2.3),
- Qty 1 USB Key with "LogViewer" software included (CODE 61898),
- Qty.4 pipes numbered: 1, 2, EXT1 and EXT2,
- Qty.2 plastic rings 1-2 (CODE 61076)

- Three different sizes of cups for the tanks of the drawer (in total 6 cups)
 - Small size (CODE 61070), diameter 28-33mm
 - Medium size (CODE 61085), diameter 32-38mm
 - Big size (CODE 61066), diameter 42-49mm


- Qty.2 metal rings (CODE 61677)

- Qty.1 allen wrench; see the chapter 3.3.1,
- Qty.1 Metallic plate (CODE 68466 to be inserted in the cavity, see the chapter 2.7),

- Qty.1 Rapid connector with hose suitable for Politainer tanks (CODE 68630) (if necessary) – refer to Appendix F.

- Qty.10 rolls of papers for the printer (if necessary) – refer to Appendix D.1.1 for information on storage (CODE 68331/A for 1 roll),
- Qty.1 plastic stylus to work with the touch-screen terminal,
- Qty.1 Operator manual in electronic format
- Qty.10 Bag sealer with handle (CODE 68559, refer to Appendix E) (if necessary)
- Qty.1 Operator manual in paper-copy.

2.2.1. Accessories

All accessories are not included in the standard delivery set and must be ordered separately.

1. External exhaust FORMALIN fumes filter kit (CODE 61550/F)

External exhaust FORMALIN fumes filter kit includes:

- Metal housing - 1 pc.,
- Filter - 1 pc.,
- Labels to the filter - 1 pc.,
- Kit of 3 different lengths exhaust tubes complete with metal clamps - 1 pc.

See detail information in chapter 4.1.

2. External Self Standing XL exhaust FORMALIN fumes filter kit (CODE 68620);

External Self Standing XL exhaust FORMALIN fumes filter kit includes:

- Metal housing - 1 pc.,
- XL exhaust FORMALIN fumes filter - 2 pcs,
- Labels to the filter - 1 pc.,
- Filter dividing chamber- 1 pc.,
- Exhaust tube complete with metal clamps - 1 pc.


See detail information in chapter 4.1.

3. HEPA filter (CODE 68290); see chapter 4.3.


2.3. How to install the bag holder

At the installation time, it's required to install the bags holder on the lateral side of the unit, on the right or on the left, as you prefer. Below is shown the procedure for the right side of the unit, repeat the same if you want to install it on the left side.


Unscrew the two black knobs from the lateral panel, shown in the image below.

Position the bag holder and screw in the two knobs.


Finally position the bags to be used first.


2.4. How to install the barcode reader

At the installation time, it's required to install the barcode reader (fixed on its support) on the lateral side of the unit, on the right or on the left, as you prefer. Below is shown the procedure for the right side of the unit, repeat the same if you want to install it on the left side.

Unscrew the two black knobs from the lateral panel, shown in the image below.

Position the barcode reader and its support and screw in the two knobs.

Finally connect the barcode reader cable to the USB port of SealSAFE.

2.5. Fume extraction system

SealSAFE is vented through the external exhaust outlet positioned on the back of the instrument, with a hose dedicated to the external fume extraction system (if any).

The Milestone SealSAFE system dispenses with formalin and it must be connected to the fumes extraction system of the Laboratory to guarantee a healthy environment. The external fume extraction system must be properly installed at all times, according to local regulations in force on laboratory safety.

The customer is required to adapt the hose connector. The plastic hose supplied has these characteristics:

- Length: 3m $\pm$ 5%
- External diameter: 66 mm $\pm$ 5%
- Internal diameter: 63 mm $\pm$ 5%
- Weight: 340 g per meter $\pm$ 10%

The external fume extraction system must permit an air flux of at least 50 m³/h (cube meter per hour).

The unit is equipped with a HEPA filter, not accessible by the user. This type of filter is used to filter the vapors coming from the biological specimens.


The filter should be replaced every 6 months during the Preventive Maintenance action performed by authorized and trained service personnel.

The management of the HEPA filter is also performed by the software, as explained into the chapter 3.4.3.


It is possible to connect an external exhaust fume filter kit. Two types of filters are available. See chapters 4.1 and 4.2 for more information.

Insert the exhaust tube (CODE 61377) onto the proper outlet (back panel).

Use the metal tube clamp to fix the hose to the exhaust fume pipeline.

Now connect the other side of tube to the fume extractor of the laboratory.

If a fume filter kit is installed on the unit, the exhaust tube should be connected to the filter, as described in the following paragraphs. Refer to the chapter 4.1 for filter type information and on how to replace it.

2.5.1. How to install the External exhaust FORMALIN fumes filter kit


External exhaust FORMALIN fumes filter kit includes:

- Metal housing - 1 pc.,
- Filter - 1 pc.,
- Labels to the filter - 1 pc.,
- Kit of 3 different lengths exhaust tubes complete with metal clamps - 1 pc.


The CODE of the Filter kit is 61550/F. This kit is an accessory, as explained in the chapter 4.1.

Follow the procedure below to correctly mount the filter.


- Switch off the unit.
- Remove the pipe connector from the rear side of the SealSAFE unscrewing the four screw that fix it.


- Take the filter that you received and connect it with already installed connector and verify that two metal claps are fixed well.


- Fix the bracket of the filter as shown in the picture below in vertical way.

- Fix it in vertical way as shown the picture below.


Connect to the
external fumes
extraction system

DO NOT BLOCK THE FILTER OUTPUT.

The silicon stopper that you find in the kit is **NOT USED** in SealSAFE configuration.

Now connect the other side of tube to the fume extractor of the laboratory, as shown in the previous picture.

The SealSAFE unit can work without being connected to a fume extraction system, only if **ALL** the following conditions are fulfilled:

- The SealSAFE must be equipped with the External exhaust fumes filter kit CODE 61550/F.
- The filter present inside CODE 66492/F must be replaced every 3 months or 1200 cycles of use ** (equivalent to 300L of 10% formalin), as explained into the chapter 3.4.4.
- The unit is operated and maintained properly as described in section 4.1.
- Presence of air recirculation in the environment of the installation compliant with government regulations.

**The frequency of the filter replacement is estimated on the basis of standard use of the equipment, considering an average use of 0.25L of 10% formalin per specimen. Specific environmental factors such as: the presence of other reagents, chemical elements in the laboratory of installation, a prevalent use of big specimens and spills of formalin into the cavity not immediately cleaned will require more frequent replacement of the filter cartridge.

If the unit is connected to the extraction system of the laboratory, replace the filter every 6 months. If NOT connected, replace the filter every 3 months or 1200 cycles of use, as explained above.

The use of the fumes filter as specified does not replace the requirement of authorized personnel to carry out environmental monitoring, nor can it guarantee the absence of staff exposure to chemicals or the removal of personal protective equipment, as required by government regulations of the final country where the unit is installed.

2.5.2. How to install the External Self Standing XL exhaust FORMALIN fumes filter kit

External Self Standing XL exhaust FORMALIN fumes filter kit includes:

- Metal housing - 1 pc.,
- XL exhaust FORMALIN fumes filter - 2pcs,
- Labels to the filter - 1 pc.,
- Filter dividing chamber - 1 pc.,
- Exhaust tube complete with metal clamps - 1 pc.

The CODE of the Filter kit is 68620. This kit is an accessory, as explained in the chapter 4.2.

At the first installation insert 2 cartridges and their separator into the filter following the instructions and the corresponding pictures below.


1. Unscrew the 14 screws from the housing and remove the filter's housing.
2. Insert the first cartridge.
3. Insert the separator of two cartridges.
4. Insert the second cartridge.
5. Close the housing of the carter and screw carefully 14 screws.

1

2

3

4

5


Finally, make sure that the screws on the cover are securely tightened and that there are no air leaks between the cover and the base.

Follow the procedure below to correctly mount the filter.


- Switch off the unit.
- Remove the plastic cup from the output of the filter.

- Connect the SealSAFE to the filter through the pipe provided by using the metal clamps.
- Verify that two metal clamps are fixed well and that there are no air leaks.

- Leave empty area below the filter (as shown in the following picture):

- Position the filter as you prefer by moving the wheels.

This model of filter is stand alone, it should not be connected to the fume extraction system of the laboratory.

The SealSAFE unit can work without being connected to a fume extraction system, only if the following conditions are fulfilled:

- The SealSAFE must be equipped with the External Self Standing XL exhaust FORMALIN fumes filter kit CODE 68620.
- The External Self Standing XL exhaust FORMALIN fumes filter CODE 68645 includes activated charcoal cartridge (2) and filter dividing chamber (1). Two filters present inside CODE 68645 must be replaced every year or 4000 cycles of use ** (equivalent to 1000L of 10% formalin), as explained into the chapter 3.4.4.
- The unit is operated and maintained properly as described in section 4.2.
- Presence of air recirculation in the environment of the installation compliant with government regulations.

**The frequency of the filter replacement is estimated on the basis of standard use of the equipment, considering an average use of 0.25L of 10% formalin per specimen. Specific environmental factors such as: the presence of other reagents, chemical elements in the laboratory of installation, a prevalent use of big specimens and spills of formalin into the cavity not immediately cleaned will require more frequent replacement of the filter cartridge.

The use of the fumes filter as specified does not replace the requirement of authorized personnel to carry out environmental monitoring, nor can it guarantee the absence of staff exposure to chemicals or the removal of personal protective equipment, as required by government regulations of the final country where the unit is installed.

2.6. Electrical installation of SealSAFE

2.6.1. Power supply

Power supply line needs to have a ground connection. Power line protection: 230V $\pm$ 10% version Miniature circuit breaker curve C, 10A, residual current device: 30mA.


Before connecting the plug to the laboratory's socket **CHECK** that the power supply voltage of the unit matches the power supply line of the laboratory. Damages occur when wrong power supplies the instrument.

2.6.2. Back side connection plate

Look at the back side of the unit.

Description:

See also the label of the instrument that provides some important information. Always communicate the S.N. of the unit to the Customer Support.

2.7. Final placement

SealSAFE unit can be positioned in its final place of use. When the instrument is positioned, push down the lever of the wheels to lock them.


Do not push down the lock lever of the drawer's wheel.


Before moving the instrument, please make sure that all 4 wheel levers are unlocked.

Then connect:


- the main power cable into the electrical plug that is located on the back of the unit (see the picture below),


- and the plug of main power supply to the socket.

Before turning "ON" the unit, open the drawer and check if the lever to unload the reagent from the cavity is in the "CLOSE" position, as shown below. If the lever is in the "OPEN" position, it will not be possible to create vacuum in the cavity.


Then insert the metallic plate into the cavity.

It's required to insert the metallic plate into the cavity in a correct way, as explained below.

WHAT TO DO:

Check the coupling between the screw and the corresponding slot, see the picture below.

The pictures above refer to the right side of the metallic plate; repeat the same operation in the left side to correctly position it.

Always leave the metallic plate positioned in the cavity. It is required to extract the plate before cleaning and decontaminating the unit, as explained in the chapter 6.


Now you can turn the unit "ON" as explained in the following chapter.

2.8. Safety Switch ON and OFF

The main switch is located on the right lateral side of the control panel, as shown in the figure below. The same switch must be used to switch the unit off.

Before turning ON the unit, make sure that you have followed all the instructions reported above.

It is now possible to switch the unit on by pressing the main switch, as shown below.

In order to safely switch the device off, move the main switch on the lateral side of the device to the position O.

- Device is switched off.
- Disconnect from power supply.

- Disconnect accessories, if necessary.

2.9. Reagents

Complete the set up with reagents, by providing at least one tank of fixative.

It's possible to connect the unit with 4 tanks:

- 2 tanks in the drawer with a capacity up to 20L maximum,
- 2 external tanks/containers of the desired capacity.

Do not use tanks with a capacity greater than 20L into the drawer.

The external tanks should be at the same height level of the unit.

Connect the unit to tanks containing the same type of fixative. Never mix fixatives of different types.

Refer to the chapter 3.3 to load SealSAFE with reagents.

3. OPERATING SEALS SAFE


3.1. Log IN and Log OUT

As a default, when the system is switched on, the LOG IN screen appears.

Press Administrator button, insert the password (issued as a separate document with the Operator Manual) and press Enter (see figure below).

The Administrator HOME PAGE appears. As default, it shows the main functions (see figure below). In whatever screen you are, press the Home button to return to the HOME PAGE (see the circle).


You can customize the HOME PAGE under the USERS settings, choosing between Favorites and Main functions (see the following paragraph).


SERVICE Log-in is available for authorized personnel only (ADMINISTRATOR can access to a restricted SERVICE functions. SERVICE has complete access to the unit functionalities).

3.1.1. Create a new user

Press the Advanced Functions icon from the status bar. Then go into the Settings screen and select Users to create a NEW USER.


The USERS settings screen appears, as shown below. Click on "Add", a keyboard appears. Enter the NEW USER name and press Enter. In the example below the name is "lab".

Now select and configure the parameters of the NEW USER choosing:

- **USER CATEGORY:** select "User" or "Admin".
- **USER PASSWORD:** click on the password box, a number pad appears as shown below, enter the password (max 6 numbers) and confirm.

Select the appropriate FUNCTIONS for the new user:

User configuration :

☐ Set 'Favourites' as Home Page

☐ Customize processings


☐ Replace the reagents

☐ Modify reagent settings

Selecting:

"Set Favorites as Home page"	The new user has the Favorites screen as Home page.
"Customize processing":	The user is allowed to create and modify programs.
"Replace the reagents":	The user is allowed to replace the empty tanks.
"Modify reagent settings":	The user is allowed to set the tanks status (choosing between empty/fill the tank), open the Reagents list and modify the reagents details (Ex: name and bottle color).

Pressing the arrow (see figure below, circle) it is possible to choose the appropriate ACCESS LEVEL in the program startup. The following window appears:


Selecting:

Ask password for processing procedure	The software recognizes (through a password system) who is performing some operations during the processing (who starts the program, aborts the program).
Ask password for reagent management procedure	The software recognizes (through a password system) who is performing some operations during the reagent management (who replaces the empty tanks, or moves them through the Moving function).

3.2. Set Date and Time

From the Advanced Settings screen, enter into Settings and select Date/time.

In the time setting window, it is possible to set the TIME ZONE, thus activating the AUTOMATIC update of time during the year.

After changing this setting, reboot the system.

3.3. How to load SealSAFE with reagents

The SealSAFE reagent drawer hosts 2 tanks. It's possible to use commercial tanks up to 20 liters each and dimensions (LxHxW) not more than (73,5 x 50,5 x 42) cm.

Two additional reagent fill ports are located on the back side of the unit, these allow for filling from higher capacity containers ($\geq 20L$).


Before loading reagents, refer to the chapters 1.9 and 2.9.

3.3.1. Reagents set up

It is necessary to first set up the unit with reagents before operating the unit.

Select the Reagents icon on the STATUS BAR (see the arrow below).


Refer to the table below for reagents icon status:

ICON	COLOR	MEANING
	Red	All tanks are empty
	Yellow	Only 1 tank is full
	Green	More than one tank is full

At the installation time the icon is red before the setup of reagents.

When you press the reagents icon, the following screen appears.


Tanks "1" and "2" refer to the tanks to be positioned into the drawer, "EXT1" and "EXT2" refer to the external connections on the back of the device which enable for filling from large containers (drums or cubes).

To operate the unit, it's required to setup as FULL at least one tank.

For example select the tank number 1 and press "Fill the tank" (circle in figure above). A message advising that the loading date will be reset to "today" appears; press OK.


The following screen appears. Follow the instructions on the screen.


Open the drawer.

There are numbered stickers 1-2 on each pipe position. Each pipe is provided with a RAPID CONNECTOR to simplify the insertion/extraction. Just press the rapid connector to easily insert/extract any pipes from each tank, as shown below.

Start the set up from the hose number one, as in the software the tank to be inserted is the number 1.


Press the rapid connector corresponding to the tube number 1, as shown below.

Choose the plastic ring suitable for the reagent to be used in a defined position (in this case number 1) and insert it in the metallic insert.


Then insert the cap in the metallic insert. There are three types of silicon caps, suitable for tanks holes of 42- 49mm, 32-38mm and 28-33mm diameter respectively.

At the end, insert the METAL RING in the metallic insert. Screw it with the allen wrench (provided with the unit) at approximately 10cm from the bottom part of the metallic insert.

The image below shows the final configuration.


The plastic and metal rings and the caps are provided only for the 2 tanks to be inserted in the drawer.

Insert the tank containing fresh reagent in the correct position (hose number 1 in this case), then insert the rapid connector in the corresponding position (number 1).


Always check that you have properly attached the rapid connectors.

In the software, press "Done" twice, as shown below.

Before pressing "Done" for the second time, it is possible to insert the reagent lot code, as shown in the previous picture. This function is not compulsory.


If inserted, it will be possible to see the lot code in the Events log, as shown in the example below.


The reagents page is automatically opened and now tank number 1 is FULL.

You can repeat the same process for the remaining tanks, if available.

The reagents icon becomes "yellow" because there is one tank full. The loading date of the tank is saved. The loading date is important because the unit will load reagent always from the oldest tank during a "Vacuum+fixative" program, as explained in the chapter 3.7.

When more than one tank is setup as full, the reagents icon becomes "green". For example below the icon is "green" because tanks number 1 and 2 are FULL.

It's required to insert the same type of reagent in all positions.

By default, the reagent to be inserted is Formalin and its color is white.

If you want to change the name, short name and/or the color of the reagent, press the icon "Reagent", as shown below (left figure). The screen shown on the right appears. It is possible to modify the corresponding fields "Name", "Short name" or "Color" by pressing the corresponding field.


3.3.2. How to replace an empty tank

When a tank is empty, follow the instructions reported in the previous paragraph 3.3.1 to set as FULL a new tank.

3.4. Owner Settings


Enter into the Advanced functions by pressing the corresponding icon from the status bar, then select Settings and Owner. This function is available only through the Service and Administrator accounts.

The OWNER SETTINGS screen appears as shown in the picture below.

3.4.1. General

Select "General" to modify the Institute and Department name where the unit is installed.

3.4.2. Ventilations


Select "Ventilation" to manage the upper (cavity) and lower (drawer) exhaust fume fans, choosing between:

- Always ACTIVE. Both exhaust fans are activated at the switching on of the unit (suggested).
- Always OFF. Both exhaust fans are always off, even during the processing (not suggested).

OFF mode is NOT SUGGESTED.

If the unit is loaded with Formalin, the ventilation status must be "Always ACTIVE".


3.4.3. HEPA filter

Select "HEPA Filter" to manage the expiration of the filter.

Below the default limits activated:

- number of TOTAL (VO+VF) processes to replace HEPA filter: 2000
- number of MONTHS to replace HEPA filter: 6 months


It is strongly suggested **DON'T** remove the flags about HEPA filter expiration.

The button "Replace HEPA filter" should be pressed once at the installation time, to set to zero the details about expiration DATE and every time it is replaced.


When the filter is expired and the operator presses the dedicated button to start a program (VO or VF), a message appears, warning the user that it is required to contact the Service Assistance. This message doesn't block the unit and so it is possible to run the program selected.

HEPA filter needs to be replaced. Please contact your Service Assistance.

If you want to activate also a **BLOCKING MESSAGE**, the Administrator can enable the checkboxes that by default are **NOT** activated (see below).

If activated, the following blocking message appears when the operator presses the dedicated button to start a program (VO or VF).


It is not possible to use the unit because HEPA filter needs to be replaced. Please contact your Service Assistance.

When the user presses OK, the Main menu is displayed and it is not possible to use the unit until the filter is replaced.

3.4.4. FORMALIN filter


Select "FORMALIN Filter" to manage the expiration of the filter.

Below the default limits activated:

- number of TOTAL VF (Vacuum+Fixative) processes to replace the filter: 1200
- number of MONTHS to replace the filter: 3 months

Follow the expiration values reported above if you are using the External exhaust FORMALIN fumes filter kit CODE 61550/F and this filter is STAND ALONE (not connected to the fumes extraction system of the laboratory).

If the External exhaust FORMALIN fumes filter kit CODE 61550/F is connected to the fumes extraction system of the laboratory, insert the following values at the installation time:

- Nr. of MONTHS to replace FORMALIN filter: 6 months and remove all the other flags.

If you are using the External Self Standing XL exhaust FORMALIN fumes filter kit CODE 68620,

insert the following values at the installation time:


- Nr. of VF processes to replace FORMALIN filter: 4000
- Nr. of MONTHS to replace FORMALIN filter: 12 months

If the Unit is NOT equipped by a FORMALIN filter and it is directly connected to an external aspiration system, remove the flags to control the expiration of the filter.

The button "Replace FORMALIN filter" should be pressed once at the installation time, to set to zero the details about expiration DATE and every time the filter is replaced.

When the filter is expired and the operator presses the dedicated button to start a program "Vacuum+Fixative", a message appears, warning the user that it is required to replace the filter. This message doesn't block the unit and so it is possible to run the program selected. This message doesn't block the unit and so it is possible to run the program selected.

FORMALIN filter needs to be replaced


If you want to activate also a **BLOCKING MESSAGE**, the Administrator can enable the checkboxes that by default are NOT activated (see below).

☒	Nr. of VF processes to replace FORMALIN filter:	1200
☐	Nr. of VF processes that stop the unit: (if FORMALIN filter not replaced)	0
☐	Kg of fixative to replace FORMALIN filter:	0
☐	Kg of fixative that stop the unit: (if FORMALIN filter not replaced)	0
☒	Nr. of MONTHS to replace FORMALIN filter:	3
☐	Nr. of MONTHS that stop the unit: (if FORMALIN filter not replaced)	0

Expiration date: 2019-09-05

If activated, the following blocking message appears when the operator presses the dedicated button to start a program "Vacuum+Fixative".

It is not possible to use the unit because FORMALIN filter needs to be replaced.

When the user presses OK, the Main menu is displayed and it is not possible to use the unit until the filter is replaced.

See the chapters 4.1.2 and 4.2.2 to replace the filters.

It is also possible to control the expiration of the filter by the quantity (Kg) of Formalin loaded, by activating a blocking and a not-blocking message. By default they are not activated.

☒	Nr. of VF processes to replace FORMALIN filter:	1200
☐	Nr. of VF processes that stop the unit: (if FORMALIN filter not replaced)	0
☐	Kg of fixative to replace FORMALIN filter:	0
☐	Kg of fixative that stop the unit: (if FORMALIN filter not replaced)	0
☒	Nr. of MONTHS to replace FORMALIN filter:	3
☐	Nr. of MONTHS that stop the unit: (if FORMALIN filter not replaced)	0

Expiration date: 2019-09-05


3.4.5. Internet connection

Select "Internet connection" to edit and modify the network settings. When you press a numerical field, a number pad appears; enter the data and press Enter (see the image below).

Data required must be supplied by the Laboratory Network Administrator (Primary WINS not required).

For example:


IP Address	192.168.1.200
Subnet Mask	255.255.255.0
Default Gateway add.	192.168.1.1
DNS server address	192.168.1.6

3.4.6. Enables

Select "Enables" to activate:

- the printer, through the button Enable Printer (by default it is already enabled),
- the additional cavity cleaning at the end of a program VO or VF.

When the checkbox "Enable Printer" is selected, at the end of every process a label containing all the details of the process is printed (for further information please refer to the chapters 3.6 for the programs "Vacuum Only" and 3.7 and for the programs "Vacuum+Fixative").


If you don't want to use the print function and/or if the paper is not inserted in the printer, it's required to uncheck this box because the software memorizes all the labels not printed that will be printed the first time that the paper will be inserted again. Please refer to the Appendix D.1 for further information.

If the Printer is NOT active:

- the button Print last run disappears,

- the button Print daily report under "Log Book" disappears, that allows you to print the list of all the case IDs of the programs executed during the selected day.

The additional cavity cleaning by default is disabled.

When enabled, at the end of every program (VO or VF depending on which checkbox has been selected in "Owner Settings – Enables") after the sealing of the bag, an extra vacuum cycle is automatically executed to clean the air present in the cavity; then it will be possible to open the cover and to remove the bag.

This process increases the duration of every single program of about 30 seconds.


The additional cavity cleaning, even if enabled, is not automatically executed if:

- during a program the button Abort is pressed,
- the vacuum level set in the program is higher than 250mbar, to avoid that the bag inflates.

If the SealSAFE is installed in the Operating Room, it is suggested to activate the procedure (if necessary) only for VO programs, if the Unit is installed in the Anatomical Pathology, it is suggested to activate it for both programs, VO and VF.


3.4.7. Milewatch

If "Milewatch" is in use in the laboratory, an example of connection is presented below:


3.4.8. Print

Select "Print" to choose the details to be printed on a label at the end of a program. By default all is enabled.

3.4.9. Case ID

Select "Case ID" to enable all the settings related to the Case ID, as:

A. "Enable Case ID"

Both checkboxes, one for the VO "Vacuum Only" and one for the VF "Vacuum+Fixatives" programs, allow you to enable or disable the request to enter a case identifier before running a program. By default they are both enabled.

ENABLE OR DISABLE CASE ID CHECK ACCORDING TO YOUR INTERNAL LABORATORY PRACTICE (GOOD LABORATORY PRACTICE, GLP).

B. "Split Case ID"

In some cases the Case ID is composed of multiple fields of data (concatenated) and the individual fields are usually separated with a special character like ";".

It is possible to use only one of these fields, if required.

By default the checkbox is not selected.

If enabled, it is required to set:


- 1) Separator character (it is the character that splits groups of information in the Case ID),
- 2) Field of interest (it is the number of the field that will be considered into the Case ID, it is possible to choose only one field).

Below an example:

Initial Case ID: 1111; 2222


- 1) Separator character: ;
- 2) Field of interest: 1

Therefore the Case ID becomes: 1111

3.5. Run a program

After the user log in, the following screen appears.


It's possible to run 2 different types of programs by selecting the corresponding icon in the screen:

- Vacuum Only
- Vacuum + Fixative

When selecting Vacuum Only, a pre-set vacuum is applied to the vacuum bag where the specimen has been inserted (the vacuum level depends on the type of specimen). For details see the chapter 3.6.

When selecting Vacuum + Fixative, a defined quantity of fixative is loaded into the vacuum bag depending on the chosen ratio and then a pre-set vacuum is applied. For details see the chapter 3.7.

In addition you can press the icon Favorites where the favorite programs can be saved (refer to the chapter 3.11 to create Favorite programs). From this page it's possible to directly run a program by selecting the corresponding icon.

3.6. Vacuum Only program

A program "Vacuum Only" can be used in 2 different cases:

1. At the end of a surgical operation, to put under vacuum a fresh specimen before being sent to the Anatomic Pathology.
2. In Anatomic Pathology, to put under vacuum a specimen to be preserved in stock (this specimen must be well fixed).

When running this program after an organ excision, remove the excess of blood (or any other organic fluid) with a medical cloth. This avoids any putrefactive effect caused by the blood/organic fluid stagnation.

There are two ways of operation to correctly insert a specimen into a bag:

1. If two operators are available: one can keep the bag opened and the second one can insert the specimen without making dirty the external and internal upper part of a bag.
2. If only one operator is available: before inserting a specimen into a bag, reverse the upper part (as shown in the picture below) in a way to do not make dirty the external and internal upper part of a bag.

Follow the instructions below to run a "Vacuum Only" program.

Before running the program, it is strongly required to clean the upper internal of sealing surfaces of the bag using dry paper, if dirty (blood, fat, tissue components, etc.). Only the upper part should be cleaned (see figure 1 and 2).

Fig.1


Fig.2

- Press the Icon Vacuum only from the main menu and a keyboard named "BAG ID" appears, as shown below.

- The barcode data is automatically entered into the system after being scanned, as shown below. The code of the bag is positioned on its lower side, as shown in the figure below. Below the label containing the barcode of the bag is shown in detail: on the top there is the bag type and below the barcode.

When the barcode is correctly read, the reader emits a confirmation sound.

- If the function "enable Case ID" is active (refer to the chapter 3.4), when the bag ID (barcode) is read, automatically another keyboard appears, named "Case ID". Scan the ID of the patient as shown below, or enter it manually by the keyboard and press "Enter".

The maximum number of characters allowed for CASE IDs is 70.

The CASE ID can contain the following characters:


- numbers: 0123456789
- uppercase letters: ABCDEFGHIJKLMNOPQRSTUVWXYZ
- lowercase letters: abcdefghijklmnopqrstuvwxyz
- symbols: !#\$%&'()*+,-.;=@[]^_`{}~
- space.

It is not possible to use characters other than those indicated.

The SealSAFE keyboard does not contain: ' + ; = [] ` { } ~

If you enter the barcode of the bag and/or of the patient by the barcode reader (not manually) you don't have to press any button on the screen but the code is automatically read and included in the details of the program.

The following screen appears, showing the Bag ID previously scanned, the Case ID (if present) and the default vacuum protocols P1-P4.

Additional information shown in the screen is the weight measured in the cavity; it should be zero because the bag has not been inserted yet.

When the Bag ID is scanned, automatically a TARE of the weight measured into the cavity is performed.

In case it is necessary to perform the TARE again, it is possible to press the button TARE before inserting the bag into the cavity and then choosing the program of interest.


- Insert the bag into the cavity. The vacuum bag with the specimen must be placed with the open side directed toward the operator.


The plastic bag must be placed with the metallic sealing bar between the two white lines (see picture below).

Always position the topmost lines on the bag over the sealing bar.


The base of the bag should rest on the base of the metallic scale; if not, it's required to use the movable metallic tool to keep the base of the bag in position and so as to prevent it from slipping downwards. See the image below. When you don't use it leave it flat on the metallic slope.

Metallic tool is UP to support the base of the bag.


- When the bag has been correctly positioned (with the metallic sealing bar between the two white lines), lower the two magnetic sidebars that help to keep the bag correctly positioned.

- Above it's shown the procedure to lower the sidebar on the left side of the cavity, repeat the same on the right one.


In the picture below you can see an example of a bag correctly positioned.

It's very important to lower the magnetic sidebars before running a process to keep the bag in position and seal it correctly.

- Now you can read the weight of the specimen before starting the vacuum program. For example below is shown a weight of 180g.

- Choose the proper vacuum program from P1 to P4.

The right program selection must be done according to the type of specimen.

Under the icons P1-P4 there is the button "i" (as indicated by the arrow below); if pressed, the list of all the types of specimens to be used for each vacuum program is shown.


Press the icon "i" and follow the indications to choose the vacuum program according to the type of specimen that you have **ONLY** if it is fresh. If the specimen is already fixed and you want to put it under vacuum for archive, it is strongly suggested to use the program P3 for **ALL** types of specimens.

When the specimens are **FRESH**, you can also refer to the Vacuum Only Selection Guidelines to choose the proper vacuum program (Appendix B). It contains the same information that you get by pressing the info icon. If you want to run a "Vacuum Only" program of a bag containing a specimen immersed in fixative, do not use the P1 program. If you want to run a "Vacuum Only" program to put under vacuum a specimen for archive, do not use the P1 program. It is strongly suggested to use the program P3 for **ALL** types of specimens, as explained above.

The programs P1-P4 differ only in vacuum level; P1 applies the highest level of vacuum and P4 the lowest.

- When you select the vacuum program (for example P1) the following screen appears.

In the lower part of the screen appears the message "Close cover and keep pressed for few seconds" that indicates the operator to close the lid and hold it for a few seconds to facilitate the creation of vacuum in the chamber.

- Now the operator has to manually close the cover and to keep it pressed for few seconds, as shown below. When the cover is closed, the vacuum pump will start creating vacuum inside the chamber.

It's very important to keep the cover pressed for few seconds when the vacuum pump starts (figure above).


The vacuum progression in the chamber is shown on the control panel display under "Sealing time left".


While the vacuum process is running, if you press Abort (see the arrow in figure above), the program will be immediately stopped and the cover will be able to be opened.

While the vacuum process is running, if you press Stop & Seal (see the arrow in figure above), the program will be immediately stopped, the plastic bag will be sealed and the cover will be able to be opened.

- When the desired vacuum level is achieved, the vacuum bag is automatically sealed.
- On completion of the sealing procedure, an audible notice sounds and a message advises the user that is possible to open the cover.

- Automatically a label showing the process data is printed (if the printer is activated, refer to the chapter 3.4).

If a CASE ID has a number of characters >14, it will be automatically printed vertically, otherwise it will be printed horizontally.

It is possible to choose the details to be printed, as explained in the chapter 3.4.


The data printed are (if enabled):

- > PROCESSING (Name of the program)
 - > Case ID
 - > Bag ID
 - > Category ("VACUUM ONLY")
 - > Specimen weight (g)
 - > User
 - > Date
 - > Time
 - > Outcome
- Open the cover.


BEFORE EXTRACTING THE BAG FROM THE CAVITY, IT IS REQUIRED TO VERIFY THAT:

- The vacuum bag should not have moved.

- Visually check the seal: it must be homogeneous and straight. The two sealing lines must be parallel and complete.


In this picture the details of good sealing lines are highlighted with yellow arrows.


- Before removing the bag from the cavity, gently force along the sealing length to check if correctly sealed. This is required to verify that the sealing is uniform all along the sealing surface of the bag, by applying a slight tension by hand over the seal.

Wait almost 15 seconds after the sealing before testing the sealing quality (see figures A, B and C).

At the end of a sealing program, it is strongly suggested to remove the Milestone vacuum bags keeping the sealed part in the raised vertical position and not to invert the bag.


CORRECTIVE ACTION WHEN SEALING A VACUUM BAG FAILS:


1. Carefully remove the bag from the cavity of the unit, keeping the partial or unsealed part of the bag in the raised vertical position (see figure below).

2. Fully reopen the bag with scissors.
 3. To now reseal the open bag clean the upper internal surfaces of the bag, as described above.
 4. Run again the SAME vacuum program.
- Tear off the label using the jagged edge.

- Remove the protective film on the back of the label, as shown below. There is a cutting line in the centre of the film that should facilitate its removal, so it is suggested to remove first one part of the film, then the second one.

- Attach the label to the bag.


Do not attach the label on the top edge of the bag because it may need to be opened and sealed again. In total it's possible to open the bag and seal it up to 4 times (maximum). Please refer on Appendix D to have more information regarding the printer.


It is required to leave bags and labels away from direct sunlight.

When working with Ethanol based fixatives, it is required to not attach the label to the bag because it may discolor if in contact with that kind of fixative. Please in this case don't remove the protective film and insert the label into the plastic pocket of the bag.

If you need to print the label more times, you can press the icon "Print the last process" on the main menu.

- Now the operator can annotate the patient data, by writing on the label of the bag (if desired).
- When running this program at the end of a surgical operation, the operator requesting the histological examination will fill the "Requisition Form". This document needs to be placed in the external pocket on the left of the label where the first annotation is written. At this point, the operator can seal the requisition form in order to avoid any accidental separation from the specimen.

When running this program at the end of a surgical operation, the plastic bag containing a fresh specimen needs to be placed into an appropriate refrigerator under controlled temperature conditions (at 4°C).

3.7. Vacuum + Fixative program

A program Vacuum + Fixative can be used in 2 different cases:

1. At the end of a surgical operation, to put under vacuum a fresh specimen with a desired quantity of fixative before being sent to Anatomic Pathology (when they don't want to receive fresh specimens),
2. In the Anatomic Pathology, to put under vacuum a specimen with a desired quantity of fixative to run, for example, the overnight fixation into the bag.

Please Refer to the Appendix A for the Vacuum+Fixative guidelines.

Follow the instructions below to run a Vacuum+Fixative program.

Before running the program, it is strongly required to clean the upper internal of sealing surfaces of the bag using dry paper, if dirty (blood, fat, tissue components, etc.). Only the upper part should be cleaned (see figure 1 and 2).

Fig.1

Fig.2

- Press the icon Vacuum+Fixative from the main menu and a keyboard named "BAG ID" appears, as shown below.

- The barcode data is automatically entered into the system after being scanned, as shown below. The code of the bag is positioned on its lower side, as shown in the figure below. Below the label containing the barcode of the bag is shown in detail: on the top there is the bag type and below the barcode.

When the barcode is correctly read, the reader emits a confirmation sound.

- If the function "Enable Case ID" is active (refer to the chapter 3.4), when the bag ID is read, automatically another keyboard appears, named "Case ID". Scan the ID of the patient as shown below, or enter it manually by the keyboard and press "Enter".

The maximum number of characters allowed for CASE IDs is 70.

The CASE ID can contain the following characters:

- numbers: 0123456789
- uppercase letters: ABCDEFGHIJKLMNOPQRSTUVWXYZ
- lowercase letters: abcdefghijklmnopqrstuvwxyz
- symbols: !#\$%&'()+,-.;=@ [^_`{}~
- space.

It is not possible to use characters other than those indicated.

The SealSAFE keyboard does not contain: ' + ; = [] ` { } ~

If you enter the barcode of the bag and/or of the patient by the barcode reader (not manually) you don't have to press any button on the screen but the code is automatically read and included in the details of the program.

- The following screen appears, showing the Bag ID previously scanned, the Case ID (if present) and the default vacuum protocols. The default ratios "weight of specimen : weight of fixative" are: "1:1", "1:2", "1:2.5" and "1.3", as shown below.

If you want to apply a different ratio, you can customize a default program under "Programs" function, as explained into the chapter 3.10.

Milestone s.r.l does not suggest any ratio. An internal protocol in the laboratory/Institute should be validated.

When the Bag ID is scanned, automatically a TARE of the weight measured into the cavity is performed.

In case it is necessary to perform the TARE again, it is possible to press the button TARE before inserting the bag into the cavity and then choosing the program of interest.

- Insert the bag into the cavity. The vacuum bag with the specimen must be placed with the open side directed toward the operator.


When running the program after an organ excision, remove the excess of blood (or any other organic fluid) with a medical cloth. This avoids any putrefactive effect caused by the blood/organic fluid stagnation.

Before running a Vacuum + Fixative program in Anatomical pathology, it's required to grossly dissect large organs and to open natural cavities for proper fixative penetration.


The plastic bag must be placed with the metallic sealing bar between the two white lines (see picture below).

Always position the topmost lines on the bag over the sealing bar.


The base of the bag should rest on the base of the metallic scale; if not, it's required to use the movable metallic tool to keep the base of the bag in position and so as to prevent it from slipping downwards. See the image below. When you don't use it leave it flat on the metallic slope.


Metallic tool is UP to support the base of the bag.

- When the bag has been correctly positioned (with the metallic sealing bar between the two white lines), lower the two magnetic sidebars that help to keep the bag correctly positioned.

- Above it's shown the procedure to lower the sidebar on the left side of the cavity, repeat the same on the right one.


In the picture below you can see an example of a bag correctly positioned.

It's very important to lower the magnetic sidebars before running a process to keep the bag in position and seal it correctly.


- Now you can read the weight of the specimen before starting the vacuum program.
- Select the Vacuum+Fixative program depending on the desired ratio "weight of specimen : weight of fixative", for example "1:2". The following screen appears.

- In the lower part of the screen appears the message "Arm injection nozzle".
- Now the operator can insert the injection nozzle into the bag by rotating 90° clockwise until it stops and locks into position.

- Check again that the open side of the bag is positioned on the sealing bar and make sure that the edge of the bag protrudes approximately 1 cm, as shown below.

- Position the upper & lower sides of the bag over the sealing bar and insert your hand inside the bag to flatten the material under the nozzle against the metallic slope. It is **very important** that the 2 sides of the bags are not attached to each other and that the lower side is flat.

- Now the operator has to manually close the cover and keep it pressed for few seconds until the software progress bar begins, this will ensure that vacuum is correctly produced in the chamber.


In the example reported in this chapter, the ratio "1:2" is selected and it means that if the weight of specimen is 460g, the unit will load into the bag 920g of fixative.

By default each type of bag has a minimum and maximum limit that controls the filling of reagent:


1. **MINIMUM LIMIT:** it is the minimum quantity of fixative that will be loaded into the bag, even if the real quantity of fixative to be loaded should be LOWER (if calculated on the base of the weight of the specimen and the selected ratio); this value avoids cases when the specimen is too small for the type of bag used.
2. **MAXIMUM LIMIT:** it is the maximum weight allowed into the bag (fixative weight + specimens weight), even if the real quantity of fixative to be loaded should be HIGHER (calculated on the base of the weight of the specimen and the selected ratio); this value avoids cases when the fixative to be loaded is too high for the volume of the bag and it may leak out from the bag.

The MINIMUM and the MAXIMUM limit values depend on the type of bag; before running a program when the user scans the BAG ID of the bag, the software will identify which is the type of bag in use and by consequence its limits. If the weight of the reagent to be loaded is less than the MINIMUM quantity established for the bag in use, a message appears advising the user that the minimum quantity will be loaded. Press "OK" to continue. In the figure below the value "40g" is only an example because the minimum quantity depends on the size of the bag.

The fixative weight declared does not include to quantity of liquid (more or less 60g) loaded into the bag at the end of a procedure, when the fluid injector is cleaned before going out of the bag.

If the sum between the weight of the fixative to be loaded (depending on the ratio chosen) and the weight of the specimen is higher than the MAXIMUM limit defined for the bag, a message warns the user that it is not possible to reach the selected ratio. The new ratio is shown. Press "Yes" to continue or "No" to abort the program. Below an example.

In the table below there are the MINIMUM and MAXIMUM limits for all types of bags:

DESCRIPTION	MINIMUM LIMIT (Fixative weight)	MAXIMUM LIMIT (Fixative weight + Specimen weight)
Vacuum bags SU LARGE	150 g	3000 g
Vacuum bags SU SMALL	40 g	300 g
Vacuum bags EXTRA LARGE	150 g	3000 g
Vacuum bags LARGE	100 g	2000 g
Vacuum bags MEDIUM	40 g	150 g
Vacuum bags SMALL	40 g	40 g

For Vacuum bags SMALL, weight of the fixative (minimum limit) and weight of the fixative + specimen weight (maximum limit) are equal, since the biological sample (biopsy) placed in this bag has a small weight - a few milligrams.

Allowed deviation range for weighing ± 10 g.

If the weight of the reagent to be loaded is higher than the MAXIMUM quantity established for the bag in use and the maximum ratio that is possible to reach is less than 1:0.2, a message appears advising the user that the specimen is too heavy and that it's not possible to continue.

In the 2 cases above, it's required to manage the specimen in a different way:


- Use a bigger bag (if available) or
- Split the specimen into 2 bags or
- Insert the specimen in an external plastic container to execute the fixation.

When the cover is closed, the vacuum pump will start creating vacuum inside the chamber.

The vacuum progression in the chamber is shown on the control panel display under "Sealing time left".

It's very important to keep the cover pressed for few seconds when the vacuum pump starts (see below).


What happens during a "Vacuum + Fixative" program, when the cavity is closed:

- A fixed level of vacuum is applied, that is necessary to load the reagent into the bag.
- After a few seconds the vacuum pump stops, some air is pumped into the bag to open it a bit and then the reagent is loaded into the bag.
- During the loading step, the "Sealing time left" bar disappears and the Fixative weight and the message "Loading reagent, do not open the cover" are shown.


- The loading stops when the desired RATIO "weight of specimen : weight of fixative" is achieved (in the example below the specimen weight is 460g, the ratio is 1:2 and so the loading stops when the fixative weight is 920g).
- For a few seconds the fluid injector sprinkles the remaining reagent into the bag.
- Then it comes back automatically to the initial position (out of the bag).
- Now the vacuum pump restarts to reach a fixed level of vacuum.

While the vacuum process is running, if you press Abort (see the arrow in figure above), the program will be immediately stopped and the cover will be able to be opened.

If you press "Stop & Seal" (see the arrow in figure above), the program will be immediately stopped, the plastic bag will be sealed and the cover will be able to be opened.

- Finally the bag is automatically sealed.
- On completion of the sealing procedure, an audible notice sounds and a message advises the user that it's possible to open the cover.
- Automatically a label showing the process data is printed (if the printer is activated).

If a CASE ID has a number of characters >14, it will be automatically printed vertically, otherwise it will be printed horizontally.

It is possible to choose the details to be printed, as explained in the chapter 3.4.


The data printed are (if enabled):

- PROCESSING
 - Case ID
 - Bag ID
 - Category ("VACUUM+FIXATIVE")
 - Specimen weight (g)
 - Fixative weight (g)
 - Spec/Fix ratio (Specimen weight / Fixative weight ratio. If the selected ratio cannot be reached because of the MIN or the MAX limits, the "Set ratio" and the "Effective ratio" will be shown)
 - Reagent
 - User
 - Date
 - Time
 - Outcome
- Open the cover.


BEFORE EXTRACTING THE BAG FROM THE CAVITY, IT IS REQUIRED TO VERIFY THAT:

- The vacuum bag should not have moved.


- Visually check the seal: it must be homogeneous and straight. The two sealing lines must be parallel and complete.


In this picture the details of good sealing lines are highlighted with yellow arrows.

- Before removing the bag from the cavity, gently force along the sealing length to check if correctly sealed. This is required to verify that the sealing is uniform all along the sealing surface of the bag, by applying a slight tension by hand over the seal. Wait almost 15 seconds after the sealing before testing the sealing quality (see figures A, B and C).

At the end of a sealing program, it is strongly suggested to remove the Milestone vacuum bags keeping the sealed part in the raised vertical position and not to invert the bag.


CORRECTIVE ACTION WHEN SEALING A VACUUM BAG FAILS:


1. Carefully remove the bag from the cavity of the unit, keeping the partial or unsealed part of the bag in the raised vertical position (see figure below).

2. Fully reopen the bag with scissors.
 3. To now reseal the open bag clean the upper internal surfaces of the bag, as described above. Program to be used: P3 program (Vacuum Only). See the chapter 3.6 to get instructions on how to run a "Vacuum Only" program.
- Tear off the label using the jagged edge.

- Remove the protective film on the back of the label, as shown below. There is a cutting line in the centre of the film that should facilitate its removal, so it is suggested to remove first one part of the film, then the second one.

- Attach the label to the bag.


Do not attach the label on the top edge of the bag because it may need to be opened and sealed again. In total it's possible to open the bag and seal it up to 4 times (maximum). Please refer on Appendix D to have more information regarding the printer.

It is required to leave bags and labels away from direct sunlight.

When working with ethanol based fixatives, it is required to not attach the label to the bag because it may discolor if in contact with that kind of fixative. Please in this case don't remove the protective film and insert the label into the plastic pocket of the bag.

If you need to print the label more times, you can press the icon "Print the last process" on the main menu.


The first time that the fixative is loaded in the vacuum bag, please select the labels of the Globally Harmonized System (GHS) pictograms that refer to chemical risks for formalin use and apply them on the Vacuum bag itself. They will be helpful to apply the Safety & Health procedures when using fixative.

You can find the labels of the pictograms inside the Vacuum bag box and each box contains 200 vacuum bags enveloped in packs of 50 pcs.

Each pack of 50 vacuum bags contains:

- A cover letter about the pictograms.
- 5 sheets of labels, each one contains the labels for 10 bags.

Do NOT apply the pictogram labels on the top edge of the bag, as this may be opened and re-sealed several times (up to 4).

Depending on the type of fixative used, apply on each bag the pictograms present on the Safety Data Sheet SDS provided by your supplier.

Apply the labels of the pictograms only the first time that the fixative is loaded in the vacuum bag.


3.7.1. Reagent management during a "Vacuum + Fixative" program

The reagent is always loaded from the oldest tank (in the reagent page when a tank is set as full, the loading date is saved, as explained in the chapter 3.3.1).


During a process, if a tank becomes empty before the completion of the loading step, the unit will try to load reagent for 30 seconds. On completion of this time, if the loading is not completed and the fixative weight doesn't increase for 30 seconds, the loading tank is set automatically as empty and:

- **if another full tank exists**, the instrument starts automatically to load reagent from this tank to complete the loading step,

- if all tanks are empty, an audio alarm and a message appear warning the user that it is possible to enter in the reagents page, set up minimum one tank as FULL and then the loading step will be automatically completed by the unit. See the message below.

Press "Yes" to enter directly to the Reagent screen. Set up as FULL at least one tank (as explained in the chapter 3.3). Then press the arrow (circle below) to exit from the reagent screen and the loading procedure will restart immediately.

3.7.2. How to use a vacuum bag more times for the same specimen

If you have a SealSAFE next to the operating room and in Anatomic Pathology, the SAME bag for the SAME specimen can be used up to the storage.

The software allows for sealing each bag up to 4 times.

Below is reported an example on how to use a bag for the all life of a specimen:


- Next to the Surgery room, the vacuum bag can be sealed the first time by a program "Vacuum only" if a SealSAFE is present; in this way the specimen arrives to the pathology lab in fresh conditions.

- Then in Anatomic Pathology, if a SealSAFE is present, the same bag can be opened and sealed again by a program "Vacuum+Fixative"; in this way the overnight fixation of the specimen can be executed (as you can see in the picture the bag filled by fixative is stand up).


- Finally the bag can be emptied of fixative before the extraction of the specimen.
- See below how to empty the fixative before opening the bag. First cut only a lateral side of the bag to empty the fixative and then open completely the bag.

Waste reagent should be disposed in compliance with the local regulations.


- Then the specimen can be extracted from the bag and some biopsies are taken to be processed.
- The remaining part of the specimen (not to be processed) can be sealed for the third time for the archive, without formalin (by a program "Vacuum Only") or with new fixative (by a program "Vacuum+Fixative") if the fixation is not sure.

In the image below the archive without formalin.

If the fixation is not sure, it's required to archive specimens with formalin.


The top edge of the bags stand up contains helpful indications for areas to seal and cut. This process can occur up to three times for sealing and twice for cutting. See below.

As shown in figure above, the cutting lines are dashed.


VACUUM BAG SU LARGE:

When the bag has been sealed and cut 2 times, before inserting the bag in the cavity to seal it for the third time, it's required to use the tool of the metal plate to keep the bag in position. See below. On the metal plate there are numbered positions etched on the metal plate: you have to set the tool at the position number 6.


It's important to use the bar of the metallic plate to keep the bag correctly positioned over the sealing bar because it's shorter after two cuts.

For the first and second sealing, do not use the bar: leave it flat on the metallic plate, as shown below.

VACUUM BAG SU SMALL:

When using the STAND UP vacuum bag small size it's always required to use the tool of the metal plate to keep the bag in position. See below.

You have to set the tool at the position number 5 before running the first sealing, at the position number 4 at the second sealing and at the position number 3 at the third sealing.

When using other types of bags (the not stand up), it's required to check their size before running a program and to evaluate if the metallic tool is necessary.

3.8. Power failure during a program

In case of power failure, both in "Vacuum Only" VO and in "Vacuum+Fixative" VF programs, the SealSAFE automatically switches off and the vacuum stays in the cavity.

When the power supply comes back on after a power failure, the system automatically switches on again, the remaining vacuum is released from the cavity and the cover can be opened.

If the power failure happened before the sealing of the bag, please follow the instructions below:


- "Vacuum Only" program: in case this type of program was performed during the power failure, repeat this program again when the power comes back on.
- "Vacuum+Fixative" program: in case this type of program was performed during the power failure and some formalin has been loaded into the bag, when the power comes back on, empty the liquid in a dedicated place and repeat the same program.

3.9. Retrieving saved processing results

Enter into the "Advanced functions" by pressing the corresponding icon from the status bar, then select Log Book, as shown below.


Choose Saved runs. The following screen appears showing all the SAVED PROCESSING FILES on the left, sorted by DATE. The dates are automatically created by the software.

Select one date (in the example below "19-Nov-2015") and the process(es) saved that day will appear on the right. File name is determined by the time when the program finishes and it is followed by "VO" if the process category is "Vacuum Only", with "VF" if it is a "Vacuum + Fixative".

Click on the file name to see the details and press "Open" to visualize the graphs of the highlighted file (figure below, yellow arrow).

Press the arrow to go back to the Saved processing screen.

The DATA reported for a "Vacuum Only" program are (as shown in the right figure above):

- Run name: it's the vacuum protocol used (P1-P4)
- Category
- User
- Start date and time
- Specimen weight
- Outcome: completed with success, aborted or stopped and sealed
- Bag ID
- Case ID (if present)

When running a "Vacuum+Fixative" program, the DATA reported are the same of those shown above for a Vacuum Only program, but in addition you will get in the same page the following information:

- Fixative weight
- Specimen/Fixative Ratio
- Reagent
- Tanks

Fixative weight : 640 g
Specimen/Fixative ratio : 1:1
Reagent : Formalin
Tanks : 1

For safety, every different type of bag has a minimum and a maximum fixative weight to be loaded (as explained in the chapter 3.7); for this reason when the ratio applied is higher or lower than the set ratio (to respect the minimum and the maximum bag limits), you will get also the information of the "Effective ratio", as shown below.


Run name : 1to1
Category : VACUUM + FIXATIVE
User : Admin
Start date : Thu, 19-Nov-2015, at 13:55
Specimen weight : 50 g
Outcome : completed with success
Bag ID : 355812948405632B
Case ID : ETHI
Fixative weight : 150 g
Specimen/Fixative ratio : 1:1 (Effective: 1:3)
Reagent : Formalin
Tanks : 2

In the example reported above, the specimen weight is 50g and the desired ratio is 1:1. The bag type used requires a minimum quantity of fixative of 150g to be sure to cover the specimen, for that the fixative weight is higher than 100g and the Effective ratio applied is 1:3.

3.9.1. How to delete a saved processing

Select a Process and press "Delete" to delete ONLY the process highlighted (left figure).

Select a Date and press "Delete" to delete all the processes saved that day (right figure).


Alternatively press "Delete All" if you want to delete ALL the programs stored on SealSAFE. Press the arrow to return to the main menu.

3.9.2. Export processing data from the terminal

SealSAFE provides the documentation necessary for the report activity of the laboratory. It is possible to export the data related to the executed programs and a document that includes the events log with all the information of the case.

Insert the USB memory stick into its slots: two USB PORTS are located on the lateral sides of the Touch Control Terminal, one on the right and one on the left (see the picture below).


From "Log book" select "Events log".

This function is available only for Administrator and Service.


The following screen appears, press "Export to USB" to export the complete events log (the list of the recorded events) of the SealSAFE.


This button is available only when a USB key is inserted into the Unit.


From "Log book" select "Saved runs". The following screen appears. Select a date or a single process and press "Export" to export that date or process dates to the USB. Alternatively select "Export all" to export all the dates present in the SealSAFE.

The exported data will be automatically inserted in a folder named by the name of SealSAFE instrument followed by its serial number, as shown in the example below. In this folder there are the list of the events "LogEvents.sealsafe.s3db" and the folder "Saved" containing the executed programs.

In the folder "Saved" there are some subfolders that order the saved runs by date: in each folder named by the saving date there are the runs executed in every single day that include the saving time in their name.


If you want to look at the Events log and the data of the executed runs on the PC, please use the program LogViewer provided with the unit, following the instructions below.


3.9.3. Install and operate with the LogViewer

The USB key provided with the unit contains the LogViewer folder. Double click on the "SetupLogViewer.exe" icon to install the program.

Press "Next" in the following screens.


Then the following screen is opened. Select the checkbox "Create a desktop icon", as shown below. Then press "Next".

The checkbox "Associate s3db extension" is selected by default and it is suggest to don't modify this setting. In this way it will be also possible to open automatically the Events log by the LogViewer program, when exported.

Finally press "Install" to complete the installation.

The LogViewer software is now ready to work.

Microsoft.NET Framework 4.0 is a prerequisite to this application, the setup will install it if not present on your PC. Christal Report Runtime is a prerequisite to this application the setup will install it if not present on your PC.

Plug the USB key on the SealSAFE where all processes are stored and export them on the key, as described in chapter 3.9.2.


Open the LogViewer software (go to "Start Menu\All programs\Milestone srl\Milestone LogViewer" or click on the "LogViewer" icon created on your desktop).

Plug the USB key in your PC and "Open" a process by clicking on the relative icon (green circle).


Select the folder automatically created after the exportation of the data from SealSAFE to the USB key, as shown below.

If you want to see the saved runs, press "Saved", otherwise press "LogEvents.selasafe.s3db" to see the events log.

The "Saved" folder contains the subfolders that organize the saved runs by executing date: each subfolder contains the saved runs including the saving time in their name.

Press twice over the run of interest and the data about the run will be displayed, as shown in the example below.

Unit name: SealSAFE	
Unit name	SealSAFE
Software Version	15.11
Serial number	125199002
Processing name	P1
Category	VACUUM ONLY
User	Admin
Start date	2015/12/09, at 11:01
Specimen weight	90 g
Outcome	completed_with_success
Bag ID	387912523028372A
Bag size	SelfStanding_LARGE

Press twice on the "LogEvents.sealsafe.s3db" to open the file, as shown below.

LINE	TIME	OPERATOR	DATE	TIME	UNIT	LOG	LOG	LOG	LOG
1	11:01	Admin	2015-12-09	11:01:01	SealSAFE	Start			
2	11:01	Admin	2015-12-09	11:01:02	SealSAFE	Start			
3	11:01	Admin	2015-12-09	11:01:03	SealSAFE	Start			
4	11:01	Admin	2015-12-09	11:01:04	SealSAFE	Start			
5	11:01	Admin	2015-12-09	11:01:05	SealSAFE	Start			
6	11:01	Admin	2015-12-09	11:01:06	SealSAFE	Start			
7	11:01	Admin	2015-12-09	11:01:07	SealSAFE	Start			
8	11:01	Admin	2015-12-09	11:01:08	SealSAFE	Start			
9	11:01	Admin	2015-12-09	11:01:09	SealSAFE	Start			
10	11:01	Admin	2015-12-09	11:01:10	SealSAFE	Start			
11	11:01	Admin	2015-12-09	11:01:11	SealSAFE	Start			
12	11:01	Admin	2015-12-09	11:01:12	SealSAFE	Start			
13	11:01	Admin	2015-12-09	11:01:13	SealSAFE	Start			
14	11:01	Admin	2015-12-09	11:01:14	SealSAFE	Start			
15	11:01	Admin	2015-12-09	11:01:15	SealSAFE	Start			
16	11:01	Admin	2015-12-09	11:01:16	SealSAFE	Start			
17	11:01	Admin	2015-12-09	11:01:17	SealSAFE	Start			
18	11:01	Admin	2015-12-09	11:01:18	SealSAFE	Start			
19	11:01	Admin	2015-12-09	11:01:19	SealSAFE	Start			
20	11:01	Admin	2015-12-09	11:01:20	SealSAFE	Start			
21	11:01	Admin	2015-12-09	11:01:21	SealSAFE	Start			
22	11:01	Admin	2015-12-09	11:01:22	SealSAFE	Start			
23	11:01	Admin	2015-12-09	11:01:23	SealSAFE	Start			
24	11:01	Admin	2015-12-09	11:01:24	SealSAFE	Start			
25	11:01	Admin	2015-12-09	11:01:25	SealSAFE	Start			
26	11:01	Admin	2015-12-09	11:01:26	SealSAFE	Start			
27	11:01	Admin	2015-12-09	11:01:27	SealSAFE	Start			
28	11:01	Admin	2015-12-09	11:01:28	SealSAFE	Start			
29	11:01	Admin	2015-12-09	11:01:29	SealSAFE	Start			
30	11:01	Admin	2015-12-09	11:01:30	SealSAFE	Start			

LogViewer Version - 16.12, the safety class in accordance with the standard IEC 62304-2006: Class A: No injury or damage to health is possible.

Since the installed software is only the controlling software/ firmware of the equipment. No damage or injury is possible. There is no likelihood of causing harm by destruction. Software makes no clinical decisions.

3.10. How to modify a program


This function does NOT create a NEW PROGRAM. It is necessary to choose the default program that most fits your requirements, make a copy and modify its details.

Before making any modifications to a program, please contact our Application Department at application@milestonemedsrl.com.

Enter into the Advanced functions by pressing the corresponding icon from the status bar, then select "Programs", as shown below.

The following screen appears, listing all the existing programs: default protocols as well as protocols created by the user (if any). On the left part of the screen there are VACUUM ONLY programs, on the right there are VACUUM+FIXATIVE programs.


To modify a Default protocol, select it and press "Copy" button (figure above, arrow), re-name the program and modify the fields.

Modification of a user's protocol does NOT require copying it before making the changes required.

The "Modify" and "Delete" buttons are made available only for already copied protocols. For safety reasons, the factory pre-stored protocols cannot be directly modified or deleted.

If you "Modify" a program created by the user, it is possible to modify the vacuum level but not its name.

If you press "Copy", automatically a program named "Copy of.." is created: press on that name to modify the name by a keyboard where you can insert the new name.


PROGRAMS VACUUM ONLY

It is possible to modify the name but also, if required, to modify the vacuum level set for the default program that has been copied.

To do this operation select "Modify vacuum level" and the vacuum level of the program is shown.

Select the desired level of vacuum and press the arrow positioned below on the left to save the new program.


It is not possible to create a program having a vacuum level higher than 600mbar, otherwise it will not be possible to seal the bag. This value applies a sealing of the bag "without vacuum".

PROGRAMS VACUUM+FIXATIVE

It is possible to modify the name but, if required, it is possible to modify the set ratio of the default protocol that has been copied.

Select the new desired ratio Specimen/Fixative and press the arrow positioned below on the left to save the new program.


The checkbox "Vacuum enable" is selected by default, in this way after the loading of the fixative a vacuum phase is always applied.

It is suggested to keep selected this checkbox, otherwise after the loading of the fixative the bag will be automatically sealed leaving some air into the bag.


3.11. How to set up a Favorite program

Insert FAVORITES programs in the HOME PAGE as follows:

- Press the Advanced Functions icon from the status bar. Then go into the "Settings" screen and select "Set favorites".


- Press "Not Defined" button under "Choose the Favorite to Set".

- Under "Program mode" select "Vacuum Only" or "Vacuum + Fixative" depending on the category of program that you want to set as favorite (see the circle in figure below).

- When you select the "Program mode", a list of the existing programs appears on the right under "Program". Select the program of your interest. You can change the BACKGROUND COLOR by pressing "Choose the Button Color", (see the example below).

- Press "Save Favorite" to save changes (figure above, yellow arrow).
- The program icon is now listed under "Favorite" screen. It's possible to customize the Main Menu as Favorite screen, as explained in the chapter 3.1.1.

3.12. Temperature Monitoring System

The temperature monitoring system is an independent device, it is not a medical device and is not included in the delivery set of SealSAFE.


To operate the device it is strongly recommended to use the Temperature monitoring system manufactured by MILESTONE S.r.l., Italy.
Use of systems not supplied by Milestone may impair safe operation and may void the manufacturer's warranty.

The temperature monitoring system ensures permanent monitoring of temperature-sensitive specimens and enables fast and simple documentation of the transport process to the pathology laboratory.

Temperature Monitoring System is composed of:

- CODE 68500 Data Logger Card (qty. 10 recommended)
- CODE 68550 Data Logger card desktop reader (qty. 2 recommended) including CODE 68539 Pen drive with data logger software.

Cards
CODE 68500

Card reader
CODE 68550

Main characteristics:

- Resolution: $\pm 1^{\circ}\text{C}$
- Temperature read every 2/5/10 minutes
- Operating temperature: from -10°C to 40°C
- Non-rechargeable card battery
- Battery lifespan: 2 years
- Card dimensions: 80x40x12 mm
- Card reader dimensions (USB-cable of 1.5 m): 115x45x27 mm
- Maximum number of digits in case ID: up to 240
- Maximum number of placing into the transport: 1000
- Time to reach transport temperature: from ambient temperature ($22-25^{\circ}\text{C}$) to 4°C : 30 minutes.

If you have a SealSAFE next to the operating room and it has integrated slot for temperature monitoring system, the instrument can be used to activate the cards in order to monitor the transport temperature of the samples to Pathology.


How to use the system

The system is composed by 2 stations:

- first: to activate the cards in the operating room where the specimens to be sent are collected, and
- second: to download the recorded data from the cards in the pathology laboratory, where the specimens are received.

3.12.1. Sender Mode for SealSAFE - Sending specimens


It is possible to activate the cards directly from SealSAFE. Insert the card into the slot positioned in the lateral side of the unit as shown below.

From the main screen of SealSAFE, press the icon "Temperature logger", as shown below.

Automatically the following window appears.

In the upper part of the screen the following information are provided:


- Unit ID: Serial Number of the unit
- Status: shows the status of the card that has been inserted in the slot of SealSAFE
- Battery level: it should be "Good" to use the card.
- Firmware: it shows the firmware of the card.


When the battery level displays "Low", it is recommended that card use being discontinued. The card may not contain sufficient power to complete a full temperature recording.

3.12.2. How to enter Case IDs manually


Before activating the card, it is possible to enter the codes of the samples to be sent together with the card. To use this function, press "Modify bar codes", and then "Add code", as shown below.


The keyboard appears automatically. Enter the Case ID through a barcode reader or manually, and press Enter. Repeat this operation for all Case IDs available. There are not a limit in the number of codes that can be entered.

If you enter Case IDs through barcode reader, they are automatically included in the list (the picture below shows the Case ID in the list), and if you enter them manually, press Enter on the keyboard. Case IDs are displayed in the central part of the screen as shown below.

To delete the Case ID, select it and press Delete code. Press Clear list to delete all the Case IDs as shown by the arrow in the previous picture. To return to the main screen, press the arrow circled in the picture above.


3.12.3. How to enter Case IDs automatically


It is not necessary to keep the card inserted into the SealSAFE unit during the automatic insertion of the codes, as the codes will be saved on the dedicated screen (it is nevertheless important not to shut off the instrument in order to avoid losing the codes). For example, the card could be kept in the refrigerator and only inserted into the SealSAFE unit at the time of activation.

Make sure that the following functions are activated in window Settings / Owner / Case ID:

When starting the program, after entering the case ID, a message will appear for the user asking if he wants to enter this case ID on the Temperature Data Logger.



- When "Split case ID" is enabled, only the pertinent part of the Case ID will be taken into account. In order to do this, it is necessary to enter the desired values into the fields for the "Separator character" (character used to separate the code) and the "Field of interest" (that which will be considered with respect to the separator character). An example is provided below:

After having entered the following values:

- Separator character: ";"
- Field of interest: "1"

If the following Code is entered: 1111;2222
The Case ID will become: 1111

After sealing the samples, press the icon "Temperature logger" on the main screen, as shown below.

Automatically Case IDs are displayed.

3.12.4. How to activate a card

Before activating the card, you need to choose how often the card should read the temperature in section "Read temperature every". This value should be chosen depending on the distance between the operating room and pathology laboratory.

- 2 minutes: up to 30 minutes of temperature registration
- 5 minutes: up to 3 days of temperature registration
- 10 minutes: up to 6 days of temperature registration

After entering one or more Case IDs and selecting the temperature reading frequency, if you need to continue working with SealSAFE programs, press the arrow icon on the left side of the screen to return to the main software page. At the same time, the previously entered case IDs will be stored in memory and will be visible when you return to the temperature monitoring page. Then start all SealSAFE programs and return to the temperature monitoring page to enter the missing case IDs. During these actions, the card must be in the slot of SealSAFE.


After loading all case IDs, perform the following actions:

1. Press "**Start**" to activate the temperature recording



2. If the printer is not active, a message will appear asking the user if he would like to print a card activation report.

If so, a report will be printed and handed along with the specimens during transportation from the operating room to the pathology laboratory:

3. As soon as the card is activated, a beep will sound confirming its activation and the status of the card will change to "Running".

It is also possible to activate the card in a "Delay mode": the temperature monitoring will start later, at the date and time chosen by the user.


So all the steps explained above on how to load case IDs are the same, it will be different only the way to activate the card, as explained below.

3.12.5. How to activate a card in a Delay mode

From the Temperature monitoring main screen select "Set start date", automatically the window indicated below on the right appears. Set the starting date and the time, then press "Start".

Automatically the card status is "Delayed start", as shown below. When the starting date and time set will be reached, the card status automatically changes to "Running".

The activated card must be placed in the refrigerator with specimens to be sent to the pathology laboratory.

For better temperature monitoring, Milestone recommends keeping the card away from the cooling elements of the refrigerator.


3.12.6. Receiver Mode - Receiving specimens

Read carefully the instructions reported below on how to download the cards data in both cases: when a card has been activated directly by a SealSAFE or by a computer.


Before using the software on the computer to read the card data, it's necessary to set its operating mode as "Receiver", as explained below.

Follow the instructions below.

1. Launch the software "Temperature Data Logger" by clicking the corresponding icon on your desktop.
2. In the upper part of the screen select "Tools" - "Settings".

Enter the password to change the settings (provided in a separate document).

The following window appears.
To select the mode, check the box next to "Receiver"
(see the circle in the figure).


In the settings window, it is possible to insert the following information:

- **Institute name / Department**
- **Out directories / Third Party Integration as explained below:**


1. Change the default path by selecting it and double-clicking.

2. Select the storage folder and press **Open**, so the first path will be saved.


3. It is possible to enter up to 3 different storage paths: select the second path and double-click. Select the storage folder and press **Open**


Then press **"Save and close"** to save the settings.


Message **"The new operative mode will be applied at the next application's start up"** will appear. Press OK.

To effect all settings, exit the software by pressing **"File"** → **"Exit"**, and launch again the software by pressing the icon on your desktop.


Now the software is ready to be used.

How to check case IDs obtained with SealSAFE

Cards automatically activated with SealSAFE contain more information than cards activated with a computer in the sending area (e.g. operating room).

For each case ID, the report will contain the following information:

User: user who runs the program

BagID: bag identification number

Program: program name

Date / Time: Date / Time

Sample Weight (g) [Sample weight (g)] (ONLY for SealSAFE)

Fixative Weight (g) [Fixative weight (g)] (ONLY for SealSAFE)

4. ACCESSORIES

All accessories are not included in the delivery set and must be ordered separately. See the following paragraphs.

4.1. External exhaust FORMALIN fumes filter kit (CODE 61550/F)

External exhaust FORMALIN fumes filter kit includes:


- Metal housing - 1 pc.,
- Filter - 1 pc.,
- Labels to the filter - 1 pc.,
- Kit of 3 different lengths exhaust tubes complete with metal clamps - 1 pc.

The FORMALIN filter kit is required, if the SealSAFE is NOT connected to the fume extraction system.

Refer to the chapter 2.5.1 for the filter expiration.

Below a picture of the filter.

Fan air flow range: 31–37m³/h (cube meter per hour).

The external fume extraction system must have an air flux of at least 37m³/h (cube meter per hour).

4.1.1. Handling and storage

HANDLING

Use protective gloves. Avoid breathing dust or wear a proper mask.

STORAGE

Store indoors in a place cool and dry.

4.1.2. How to change the filter

The filter cartridge has the code 61550/F.

Record on the dedicated label the date of installation and replacement/expiration of the filter cartridge to know when it is expired. Attach the label to the filter.

The label should be attached to the filter.

The management of the FORMALIN filter is done also by the software, as explained into the chapter 3.4.4.

Refer to the previous paragraph 4.1.1 before handling the filter.

Refer to the chapter 11.5 for waste disposal procedure.


Follow the instructions and the corresponding pictures below.

1. Unscrew the six black knobs from the housing.
2. Then remove the filter's housing.
3. Now remove the black seal.
4. Remove the filter.
5. Replace the filter with the new one, place the black seal.
6. Close the housing of the carter and screw the six knobs.

Then reconnect the filter to the external fumes extraction system (if not used stand alone).

Then reset the expiration date of the filter into the software, by pressing the button "Replace FORMALIN filter", as shown in the chapter 3.4.4.

4.2. External Self Standing XL exhaust FORMALIN fumes filter kit (CODE 68620)

External Self Standing XL exhaust FORMALIN fumes filter kit includes:

- Metal housing - 1 pc.,
- XL exhaust FORMALIN fumes filter - 2pcs,
- Labels to the filter - 1 pc.,
- Filter dividing chamber - 1 pc.,
- Exhaust tube complete with metal clamps - 1 pc.

The External Self Standing XL exhaust FORMALIN fumes filter kit is required, if the SealSAFE is NOT connected to the fume extraction system. Refer to the chapter 2.5.2 for the filter expiration.

Below a picture of the filter.

When working with this filter, SealSAFE should not be connected to the fumes extraction system of the laboratory.

4.2.1. Handling and storage

HANDLING

Use protective gloves. Avoid breathing dust or wear a proper mask.


STORAGE

Store indoors in a place cool and dry.

4.2.2. How to change the filter

This type of filter includes the External Self Standing XL exhaust FORMALIN fumes filter kit CODE 68645 with activated charcoal cartridge (2) and filter dividing chamber (1).

Record on the dedicated label the date of installation and replacement/expiration of the filter cartridge to know when it is expired. Attach the label to the filter.

The label should be attached to the filter.

The management of the FORMALIN filter is done also by the software, as explained into the chapter 3.4.4.


Refer to the previous paragraph 4.2.1 before handling the filter.

Refer to the chapter 11.5 for waste disposal procedure.


Follow the instructions and the corresponding pictures below.

1. Unscrew the 14 screws from the housing and remove the filter's housing.
2. Remove the first cartridge.
3. Remove the separator between the two cartridges.
4. Remove the second cartridge.
5. Replace the two cartridges with the new ones, remember to insert the separator between them.
6. Close the housing and screw carefully the 14 screws.

In the end verify the correct closure of the screws of the cover and that there aren't air leaks between cover and housing.

Then reset the expiration date of the filter into the software, by pressing the button "Replace FORMALIN filter", as shown in the chapter 3.4.4.


4.3 HEPA filter (CODE 68290)

The HEPA filter is designed to prevent environmental pollution from microbiological contamination and it is used to filter vapors from biological samples.


The HEPA filter consists of 3 main components:

- FILTER ELEMENT (Cartridge)
- HOUSING
- DRAIN connected to the housing

The filter is mounted vertically with the drain pointing downwards.

The image below shows the filter cartridge.

HEPA filter is not accessible by the user.

The filter should be replaced by authorized service personnel every 6 months of use during the Preventive Maintenance.

The management of the HEPA filter is also performed by the software, as explained in the chapter 3.4.3.

5. CONSUMABLES

All consumables are not included in the delivery set and must be ordered separately.

All the parts described in this chapter should be provided only by the manufacturer Milestone s.r.l., Italy.

5.1. Vacuum bags

It is recommended to use the medical device "Bags for TISSUE SAFE and SealSAFE units for vacuum sealing the biological samples", registration certificate in the Russia No. P3H 2019/8538, as a consumable for working with the instrument.

The particular design and the special multi-layer of MILESTONE Vacuum bags, in combination with a dedicated sealing bar (supplied by an electrical intensity specific for thickness and material), allow a perfect double sealing, that resists even in case of accidental fall or external compression.

The use of MILESTONE vacuum bags is to be considered exclusive to the SealSAFE instrumentation. The use of such bags for a different type of instrumentation or the use of vacuum bags not provided by MILESTONE for the above-mentioned instruments preclude the safety of use, resulting in biological risk for the operator and risk of impairment of quality for the specimen, making void the manufacturer's warranty and the CE-IVD.

Each bag is marked with the symbol "Single use only" (see the circle in figure below); this is intended that the bag must be used for the same specimen because it can be sealed 4 times if necessary, as explained in the chapter 3.7.2.


Milestone S.r.l. provides 2 types of stand up vacuum bags and four types not stand up. See the following paragraphs.

5.1.1. Storage of Vacuum bags

Store in a cool, dry place, away from heat, do not kept in warm storage areas or direct sunlight.

5.1.2. Stand up vacuum bags

Milestone supplies two types of Stand up vacuum bags, they differ only in size. See the image below.

REF 68409SS

REF 68545SS


REF	DESCRIPTION
68409SS	VACUUM BAG SU LARGE
68545SS	VACUUM BAG SU SMALL

Characteristics:

- Disposable bags (if necessary prescribed the sealing of the bag up to 4 times (maximum) for the same specimen from the same patient) as described in chapter 3.7.2),
- Due to the special design with the stand up bottom, the Vacuum bag is Stand Up and can stand alone on a table or other horizontal surface,
- High mechanical strength,
- Applied sealable pocket for documents that can be opened and closed more times,
- Continuous lines that suggest where sealing,
- Dotted lines that suggest where cutting,
- Double barrier layer of material BOPA/PE (polyamide and polyethylene),
- Indication of lot number and date expiration,
- Enveloped in packs of 50 pcs.

5.1.3. Vacuum bags not Stand up

Four models are available:

REF	DESCRIPTION
68047SS	VACUUM BAG SMALL
68055SS	VACUUM BAG MEDIUM
68054SS	VACUUM BAG LARGE
68046SS	VACUUM BAG EXTRA LARGE

The four types of vacuum bags described above differ only in size. Their common characteristics are:

- Disposable bags (if necessary prescribed the sealing of the bag up to 4 times (maximum) for the same specimen from the same patient) as described in chapter 3.7.2,
- High mechanical strength,
- Applied sealable pocket for documents that can be opened and closed more times,
- Continuous lines that suggest where sealing,
- Double barrier layer of material BOPA/PE (polyamide and polyethylene),
- Indication of lot number and date expiration,
- Enveloped in packs of 50 pcs.

For more information see Instructions for Use of "Bags for TissueSAFE and SealSAFE units for vacuum sealing the biological samples".

6. ADDITIONAL DEVICES THAT ARE NOT PART OF THE SYSTEM, BUT WHICH CAN BE USED WITH IT

6.1. Kit Refrigerator box

To maintain the temperature regime during transportation of biological samples, it is recommended to use additional devices "Kit Refrigerator box" or thermal bag with cool packs (registration certificate in the RF No. P3H 2016/4290 or similar products duly registered in the territory of Russia).

Specimens should be transported strictly in accordance with the Specimens Transportation Guidelines (see Appendix C).

Fresh specimens under vacuum must be stored at a temperature of +4°C.

For optimal archiving of samples, Milestone S.r.l. recommends to keep cool packs in a refrigerator next to SealSAFE at + 4 °C ± 2 °C.

Do not place or transport cool packs in Kit Refrigerator box at temperatures lower than 2°C. Frozen packs could damage the samples during transport.

It is very important to maintain the correct temperature inside the Kit Refrigerator box. Never transport any samples without cool packs present.

When the Courier is ready to transport the specimens to the Anatomic Pathology Department, the vacuum sealed specimens must be placed in the Kit Refrigerator box with cool packs. Number of cool packs should be in accordance with the inner volume of the kit refrigerator box used and number of biological specimens transported. Box should be closed tightly and then transported at +4°C.

After transportation, specimens are ready for grossing procedure.

6.2. Vacuum bags

It is recommended to use the medical device "Bags for TISSUE SAFE and SealSAFE units for vacuum sealing the biological samples" manufactured by Milestone s.r.l., Italy, registration certificate in the RF No. P3H 2019/8538, as a consumable for working with the instrument.

See detail information in chapter 5.1.

6.3. Temperature Monitoring System

SealSAFE WDLR Automatic Fixative Filling and Vacuum Sealing System has an integrated slot for temperature monitoring system. This product is not a medical device and not included in the delivery set of the instrument.

To operate the device it is strongly recommended to use the Temperature monitoring system manufactured by MILESTONE S.r.l., Italy.

Use of systems not supplied by Milestone may impair safe operation and may void the manufacturer's warranty.

See detail information in chapter 3.12.

6.4. Compatible reagents for the fixation

Information on compatible reagents for the fixation is given in chapter 1.9.

7. MAINTENANCE

Maintenance of SealSAFE must be done at two levels.

- 1) The first is at the laboratory level: simple routine cleaning maintenance can prevent the likelihood of more serious problems occurring at a later stage.
- 2) The second is technically more detailed and is beyond the scope of normal laboratory personnel.

Attention: check with local distributor or with Milestone customer support (customersupport@milestonemedsrl.com) to perform Preventive Maintenance and regular inspections in order to maintain your SealSAFE within safe operative standards.

7.1. Operating room personnel

Minimal routine maintenance is required by Operating Room personnel. The following is recommended to ensure optimal performance of the SealSAFE.

The operator maintenance can be divided in 2 categories:

- 1) ORDINARY MAINTENANCE
- 2) EXTRAORDINARY MAINTENANCE

7.1.1. Ordinary maintenance

Daily

Decontamination and Cleaning

Follow the procedure described in the chapter 12 to decontaminate and clean the following parts of the instrument:

- 1) Metallic scale
- 2) Cavity
- 3) Chamber drain hole
- 4) Cover
- 5) External surfaces

Cleaning and decontamination is carried out with chemicals.

It is recommended to use disinfectants containing polyphenols for the decontamination.

For the Russian Federation:

In accordance with the rules established by "MU 287-113" (taking into account "GOST 177", "GOST 25644"), it is allowed to clean the working surfaces with a 3% solution of hydrogen peroxide (for example, produced by JSC "Samaramedprom", Russia) with the addition of 0,5% "Lotos" type detergent (for example, produced by "RHP" LLC, Russia). It is allowed to use cleaning and disinfecting agents from other manufacturers that are similar in composition and allowed for circulation in the territory of the Russian Federation.

DO NOT use:

- Alcohol-based detergents for the cover.
- Chlorine-based solutions (e.g. bleach) at any concentration and temperature.
- Corrosive substances (acids in general) at any concentration and temperature or abrasives and sharp instruments that may damage the coating.

Check the cover and the edge surface of the vacuum chamber.

Traces of organic residues must be removed to ensure the optimum closing of the cover.

Cleaning of the drain hole in the cavity

Please do it carefully.


Flush any residual debris in the vacuum chamber drain hole with 1 liter of warm water. Please see guidance in "Extraordinary Maintenance" section 7.1.2.

Weekly**Vacuum Pump Cleaning**


It's required to run a "Pump Cleaning" process two to three times every week. To run this cycle you don't need any bag or reagent.

The procedure takes only ten minutes and it will improve the life of the vacuum pump.

From the main menu press "Pump Cleaning" icon.

**Cleaning of the Metallic screen**

Remove the metallic scale and verify that the metallic screen covering the drain hole located at the right rear corner of the cavity floor is clean. Otherwise unscrew it, wash under a stream of water and screw it again.

**When required****Roll of paper replacement (if printer activated)**

Replace the paper roll when it is completely used, as explained in the Appendix D.1.2.

Filters replacement

Replace the Formalin filter when the instrument warns the user at the beginning of a program (see the chapters 4.1.2 and 4.2.2 depending on the type of filter installed).

The HEPA filter will be replaced by Service personnel.

Yearly


Contact your technical support for a preventative maintenance contract.

7.1.2. Extraordinary maintenance**Accidental spillage of fixative into the cavity**

Refer to the chapter 3.7 for information on how to correctly position a bag into the cavity to AVOID any spillage of fixative.

What to do in case of fixative into the cavity:


- Wear individual protective devices.
- Remove the metallic scale from the cavity. Before doing this operation, clean it with paper wet with warm water, if dirty.
- Open the reagent drawer and place the white rapid connector of the dedicated drain tube in position as shown in the right figure and the other side of the tube into a waste container of correct capacity (at least 3 liters).

- Flush any residual debris in the vacuum chamber drain hole with 1 liter of warm water.
- Turn the lever down into the "OPEN" position. Any fluids and debris can now be unloaded through the chamber drain hole located at the right rear corner of the cavity floor. The hole is covered with a removable metallic screen to collect large particles (figure below).

- When the procedure is complete, turn the lever to the "CLOSE" position. If the valve is not closed, the unit will NOT be able to create vacuum in the chamber.
- Press the rapid connector tab to remove the tube. You can leave the tube attached in the wall of the drawer, as shown in the figure below.

- Dispose the waste container in agreement with the government regulations.
- Clean the lateral sides of the cavity of any residual formalin.
- Dry all cavity surfaces with paper.
- Position some dry paper over the drain hole, inside the chamber, then run a "Pump cleaning" program. The paper will absorb any residual water left in the discharge canal.


- Replace the metallic scale into the cavity. Now the unit is ready to be used again.

Accidental spillage of fixative into the drawer

What to do in this case:

- Wear individual protective devices.
- Open the drawer, put a container of correct capacity under the unit, unscrew the drainage plug and unload the residual reagent as indicated in the figure below.

- Dispose the reagent unloaded accordingly with government regulations.
- Remove the fixative tanks (if present) and place them in a vented area.
- Put again the container just emptied under the unit and unload some warm water into the drawer to clean it. The quantity of water should be compatible with the capacity of the waste container.
- Dispose the water accordingly with government regulations.

- Dry all cavity surfaces with paper.
- Screw the drainage plug, position again the tanks dried by paper and close the drawer.

7.2. Service personnel

Please contact the local distributor or directly the Milestone customer support department (customersupport@milestonemedsrl.com) to receive more detailed information about the Preventive Maintenance and regular inspections in order to maintain your SealSAFE within safe operative standards.

7.3. Calibration of the instrument

These functions are available only to the SERVICE personnel and the ADMINISTRATOR.

The SealSAFE should only be calibrated by a specially trained technician. Calibration of the load cell and pressure sensor must be performed during the installation of the device and at each preventive maintenance of SealSAFE.

Calibration of the load cell

Login a Service User and Access at the Factory setting Page.

Milestone will provide the ADMINISTRATOR password after completing the appropriate training and receiving information about the serial number of your SealSAFE. The ADMINISTRATOR password is strictly confidential information.

Insert a definite weight on the load cell and check the measured weight by the unit. Acceptability range: $\pm 10\text{g}$ the known weight; then sign the calibration points.

Weight calibration		Pressure calibration		Counters T.		Counters P.		Reserved	
09 Feb 15 17:19 Service									
Load cell									
Actual	0	(g)	1593	(raw)					
Point 1	0	(g)	1594	(raw)	Acquire Point 1				
Point 2	5000	(g)	3584	(raw)	Acquire Point 2				

Point 1 _____ (g) _____ (raw)

Point 2 _____ (g) _____ (raw)

Calibration of the pressure sensor

Check the pressure calibration values. Acceptability range: $\pm 10\text{mBar}$. To calibrate the unit, follow the instruction reported on the Service Manual MM092.

09 Feb 15
17:20
Service

Weight calibration | Pressure calibration | Counters T. | Counters P. | Reserved

Current pressure : 979 (mBar) 531 (raw)

Point 1 : 982 (mBar) 532 (raw) Go to high press.
(high press.)

Point 2 : 7 (mBar) 203 (raw) Go to low press.
(low press.)

←

Point 1 _____ (g) _____ (raw)
Point 2 _____ (g) _____ (raw)

8. Responsibility of the Manufacturer

Milestone S.r.l. is the legal manufacturer that is responsible for the quality, efficiency and safety of the registered medical product for use, normal operation and possible negative results of using it, except for cases of negligence, or incorrect use through the fault of user, or force majeure.

9. Repair and service information

Maintenance and current repairs of Milestone devices is available anywhere in the world in one of the international service centers, as well as in your area. To arrange service or repair of your instrument, please contact your nearest office; updated contact information can be found at www.milestonemedsrl.com. For service organizations and technical support see details below.

Customer's Support

Tel.: +39 035 572735

E-mail: customersupport@milestonemedsrl.com

Authorized representative in the Russian Federation:

BioLine LLC

Legal address: 197022 Saint-Petersburg, Municipal district Aptekarsky ostrov, Professora Popova street 23 lit.E, structure 3H/5H, room 208, Russia

Postal address: 197022 Saint-Petersburg, Professora Popova street 23 lit.E, Russia

E-mail: service@bioline.ru

Tel: +7(812) 320-49-49

Fax: +7(812) 320-49-40

Hotline: 8-800-555-49-40.

This instrument is an in-vitro medical device and it is usually installed in a laboratory where specimens and other biological tissues are present. For your safety, it is therefore required to clean and disinfect before entering in contact with it.

It is required in case of return that the unit is cleaned and disinfected before sending it back to MILESTONE S.R.L. Non-disinfected devices does not accepted.

In case of any undesired reactions or events during the use of the product, please immediately contact with the Customer's Support or your local distributor.

10. WARRANTY

Warranty covers the medical devices manufactured by Milestone S.r.l., Italy.

Milestone S.r.l. guarantees that the contractual product delivered has been subjected to a comprehensive quality control procedure based on international standards, and that the product is faultless and complies with all technical specifications.

The warranty period is 13 months from the date of sale.

The guaranteed shelf life of the product is at least 6 months.

This warranty becomes invalid in case of mishandling, damage and negligence and further in case of usage of inappropriate spare parts and/or accessories or unauthorized modification of the unit.

The instrument service life is 10 years.

11. WASTE DISPOSAL

11.1. Waste disposal of Vacuum bags

The Vacuum bags supplied by Milestone s.r.l. should be treated as biohazard material after their use.

Waste material should be disposed of in an approved incinerator in compliance with all federal, provincial and local government regulations.

For the Russian Federation, it is recommended to dispose the vacuum bags in accordance with the rules implemented by "SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for maintenance of the territories of urban and rural settlements, for water bodies, potable water and potable water supply, air, soils, living quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures" for the medical waste classified as class B (epidemiologically hazardous waste)..

11.2. Waste disposal of tanks

Vapors may collect in empty containers. Treat empty containers as hazardous.

Waste material should be disposed of in an approved incinerator in compliance with all federal, provincial and local government regulations.

For the Russian Federation, it is recommended to dispose the tanks in accordance with the rules implemented by "SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for maintenance of the territories of urban and rural settlements, for water bodies, potable water and potable water supply, air, soils, living quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures" for the medical waste classified as class B (epidemiologically hazardous waste).

11.3. Waste disposal of fixatives

Handle with care and store in a cool dry space using a tightly closed container.

Waste material should be disposed of in an approved incinerator in compliance with all federal, provincial and local government regulations.

For the Russian Federation, it is recommended to dispose the used fixatives in accordance with the rules implemented by "SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for maintenance of the territories of urban and rural settlements, for water bodies, potable water and potable water supply, air, soils, living quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures" for the medical waste classified as class B (epidemiologically hazardous waste).

Some reagents can have potential health effects and cause environmental pollution if not correctly disposed.

For further information please refer to the SDS (Safety Data Sheet) provided by your reagents supplier.

11.4. Waste disposal of the equipment

Before disposal, the instrument and its components are subject to mandatory decontamination (disinfection). The disinfection method is determined by the capabilities of the organization carrying out the medical activities, and is selected when developing a scheme for handling medical waste.

The unit is subject to disposal in the following cases:

- end of service life;
- in cases endangering the life and health of the service personnel.

All packaging, including transport packaging, must be disposed of.

It is prohibited to dispose of as household waste.

Decommissioning of the equipment for disposal:

- 1) Switch off the unit. Disconnect the power supply.
- 2) Disconnect all the accessories, if available.
- 3) Clean and disinfect all parts of the equipment according to the manufacturer specification.
- 4) Dispose all electric and electronic parts of the equipment.

Complying with directive 2012/19/EC of the European Parliament and of the Council of 4 July 2012 on waste electrical and electronic equipment (WEEE), the separate collection environmental managed of equipment is mandatory.

It is necessary to return the used equipment to the distributor or to inquire about the presence of a local empowered system for collection and disposal of WEEE.

The in-observance of Directive 2012/19/EC or of the local law which acknowledge it can have potential effects on the environment and human health.

This symbol indicates separate collection for electrical and electronic equipment

In case of installation outside the European Community territory, the instrument and/or its components after disinfecting must be disposed of in accordance with the existing applicable, local regulations.

For the Russian Federation, it is recommended to dispose it in accordance with the rules implemented by "SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for maintenance of the territories of urban and rural settlements, for water bodies, potable water and potable water supply, air, soils, living quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures" for the medical waste classified as class A (non-hazardous waste), excluding electronic parts. Electronic parts must be disposed of in accordance with local specified regulations with respect to household and industrial waste (in particular, office and computer technology).

11.5. Waste disposal of filters (if present)

The filter CODE 61550/F should be treated as biohazard after use.

The filter CODE 68620 should be treated as biohazard after use.

The filter CODE 68290 should be treated as biohazard after use.

Waste material should be disposed of in compliance with all federal, provincial and local government regulations.

For the Russian Federation, it is recommended to dispose the used filters in accordance with the rules implemented by "SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for maintenance of the territories of urban and rural settlements, for water bodies, potable water and potable water supply, air, soils, living quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures" for the medical waste classified as class B (epidemiologically hazardous waste).

12. DECONTAMINATION AND CLEANING

Cleaning and disinfection must be performed daily, after completion of work.

For safety reasons, any operation without the individual protective devices must be avoided. The instructions confirm the user's responsibility to carry out a proper decontamination in case of a spillage of hazardous materials above or inside the unit.

The instructions confirm the user's responsibility to NOT use decontaminating or cleaning agents that may cause danger as result of a reaction with parts of the instrument or with contaminated material inside the unit.


12.1. Safety precautions

The personnel involved in the decontamination procedure of SealSAFE and its accessories (if present) must wear the appropriate individual protective devices, in particular:


- Single use, HI-RISK protective gloves (100% hi-risk latex or equivalent). Personal Protective Equipment that meets the standard EN 374 and the Directive 89/686/CEE.

- Antibacterial mask face for surgery room and plastic apron. Personal Protective Equipment that meets the standard EN 149 Grade FFP3D and the Directive 89/686/CEE.

- Protective eyeglasses. Personal Protective Equipment that meets the standard EN 166 and the Directive 89/686/CEE.

12.2. SealSAFE decontamination and cleaning procedure

12.2.1. Detergents

Please follow the instructions below on how to carry out a correct cleaning and decontamination procedure.

Do not remove the side panels of the unit during the cleaning operations.

For the decontamination and cleaning do not use procedures different than those explained in this chapter.

Before doing any operation disconnect the power supply. Do not use direct or high-pressure water to clean the instrument.

Cleaning and decontamination is carried out with chemicals.

It is recommended to use disinfectants containing polyphenols for the decontamination.

For the Russian Federation:

In accordance with the rules established by "MU 287-113" (taking into account "GOST 177", "GOST 25644"), it is allowed to clean the working surfaces with a 3% solution of hydrogen peroxide (for example, produced by JSC "Samaramedprom", Russia) with the addition of 0,5% "Lotos" type detergent (produced by "RHP" LLC, Russia). It is allowed to use cleaning agents from other manufacturers allowed for circulation in the territory of the Russian Federation.

DO NOT use alcohol-based detergents.

Never use chlorine-based solutions at any concentration and temperature (e.g. bleach).

Do not use corrosive substances at any concentration and temperature (acids in general).

For further information, please contact the manufacturer at: application@milestonemedsl.com.

12.2.2. Procedure

1) Metallic Plate

- Remove the metallic plate from the cavity (picture below).


- Spray with disinfectant solution following the instructions reported above in the section 12.2.1 and leave it for at least 5 minutes.
- Wash under tap water to remove the disinfectant.
- Wash with an appropriate detergent and, if necessary, remove any organic trace with a non-abrasive tool.
- Wash under tap water to remove the detergent.
- Dry with paper or with a clean cloth.


2) Drain hole in the cavity

- Unload at least 1 liter of warm water into the hole present in the cavity to clean it by any residuals.


The reagent is unloaded through the hole present into the cavity, in the right side (below the metallic net), through a dedicated pipe.

- Open the reagents drawer and place the rapid connector of the dedicated drain tube in position as shown in the figure below and the other side of the tube into a waste container.
- Then turn the lever down into "OPEN" position to unload the reagent from the cavity to an external container.

- When the draining of reagent ends, press the rapid connector (see the arrow in the figure above) to remove the tube. You can leave the tube attached in the wall of the drawer as shown below.

Do not unload reagent from the cavity without connecting the tube.

- When the procedure is ended, turn the lever to the "CLOSE" position. If you don't do that, the unit will NOT be able to create vacuum in the chamber.

3) Cavity

- Follow the procedure described in the chapter 12.2.2 to unload any residual reagent from the cavity.
- Spray the cavity with a disinfectant solution following the instructions reported above in the section 12.2.1 and leave it for at least 5 minutes.
- Clean the surface with a cloth moistened with water to remove the disinfectant.
- Use an appropriate detergent to clean the cavity and to remove any organic trace with dry cloth or paper towel.
- Clean the cavity with a cloth moistened with water to remove the detergent.
- Dry the surface with paper or with a clean cloth.

During the procedure described above for the decontamination and cleaning of the cavity, pay particular attention to the sealing bar and its connectors (indicated by arrows in the figure below).

The sealing bar can be extracted from the cavity during the cleaning of it, if necessary. Then it is very important to position the bar again correctly; to do this operation push it down into the dedicated position.

Before using again the Unit, verify that the sealing bar has been positioned correctly.

4) Cover

For cleaning, we recommend soap water or a liquid detergent. Please do not use acid or alcohol based cleaners.

- Spray the internal and external surface of the plexiglass cover with a disinfectant solution following the instructions reported above in the section 12.2.1 and leave it for at least 5 minutes.
- Clean the surface with a cloth moistened with water to remove the disinfectant.
- Use an appropriate detergent to clean the cover and to remove any organic trace with dry cloth or paper towel.
- Dry the surfaces with paper or with a clean cloth.
- Dry with paper or with a clean cloth.
- Dab with particular attention the silicon bar (indicated by the symbol A in the figure below) with a dry cloth or paper towel.
- Run a P1 "Vacuum Only" program without using any bag to facilitate the drying of the silicon bar. If necessary do it a second time and then wipe with paper.
- Dry the silicon bar with paper or with a clean cloth.

During the procedure described above for the decontamination and cleaning of the cavity cover, pay particular attention to the black rubber placed on the inner surface of the cover (indicated by the symbol B in the picture above). It is important that it's clean and dry so that it fits tightly to the steel surface of the working area during the execution of a vacuum program.

5) External surfaces

- Spray the external surfaces of the unit with a disinfectant solution following the instructions reported above in the section 12.2.1 and leave it for at least 5 minutes.
- Do not spray on the back panel of the unit.
- Clean the surface with a cloth moistened with water to remove the disinfectant.
- Use an appropriate detergent to clean the surfaces.
- Clean the cavity with a cloth moistened with water to remove the detergent.
- Dry the surface with paper or with a clean cloth.

Do not spray directly on the printer, the barcode reader or into the USB slots.

6) Final steps

- Place the metallic plate in the cavity.
- Before removing gloves, wash them with soap, water and disinfectant.
- Wash your hands.
- Remove the Personal Protective Equipment.
- Make sure the contaminated material is packed and disposed in an approved incinerator in compliance with all federal, provincial and local government regulation.

On completion of the cleaning operations, wait at least 5 minutes before connecting the power supply cable. Upon completion of these procedures, the equipment is ready for use again.

13. PRECAUTIONS IN CASE OF THE MEDICAL DEVICE MALFUNCTION, FAILURE, OR DEVIATIONS IN FUNCTIONING

Before starting operation, check that the device and the mains cable are not damaged externally.

In the event of an emergency situation or equipment failure, disconnect the instrument from the power supply or follow the error messages displayed on the control panel.

Only authorized and qualified SERVICE personnel may repair the instrument and access the instrument's internal components.

In the event of a malfunction or deviations in functioning of the SealSAFE device, contact the authorized representative of the manufacturer in your region (see chapter 9) or contact Milestone customer support service (customersupport@milestonemedsrl.com) to call a technical engineer.

ERROR MESSAGE

ERROR 1 - Serial communication fault
Try to reboot the unit.

ERROR 2 - Fan not running
Contact Milestone customer support service to call a technical engineer.

ERROR 3 - Drawer fan not running
Contact Milestone customer support service to call a technical engineer.

ERROR 11 - Vacuum for load reagent not enough
Check if the black seal is positioned correctly around the cover.
Check if the HEPA filter tap is closed well,
Check if the pump is running and all the connector are connected properly without vacuum leakage,
Check if the frontal Tap inside the drawer is closed correctly
Check if the vacuum transducer and proportional valve work correctly.

ERROR 12 - Load reagent time out
Check if the black seal is positioned correctly around the cover.
Check if the HEPA filter tap is closed well,
Check if the pump is running and all the connector are connected properly without vacuum leakage,
Check if the frontal Tap inside the drawer is closed correctly
Check if the vacuum transducer and proportional valve work correctly.
Check the level in the tank from which the unit load the reagent.
Check the fast-plug-in black tube from the tank to the unit in the drawer and the seal.
Clean the pipes

ERROR 13 - Release nozzle time out
Open the dashboard and check if the Nozzle spring work correctly.

ERROR 14 - Vacuum for sample not enough
Check if the black seal is positioned correctly around the cover.
Check if the HEPA filter tap is closed well,
Check if the pump is running and all the connector are connected properly without vacuum leakage,
Check if the frontal Tap inside the drawer is closed correctly
Check if the vacuum transducer and proportional valve work correctly.

A. APPENDIX A – “VACUUM+FIXATIVE” PROGRAMS GUIDELINES

In “Vacuum+Fixative” programs, it is very important:

1. to position correctly the bag
2. to check the sealing quality at the end of a process.


A.1. How to correctly position a bag into the cavity

The vacuum bag with the specimen must be placed with the open side directed toward the operator.


The plastic bag must be placed with the metallic sealing bar between the two white lines (see picture below).

Always position the topmost lines on the bag over the sealing bar.


The base of the bag should rest on the base of the metallic scale; if not, it's required to use the movable metallic tool to keep the base of the bag in position and so as to prevent it from slipping downwards. See the image below. When you don't use it leave it flat on the metallic slope.

Metallic tool is UP to support the base of the bag.

When the bag has been correctly positioned (with the metallic sealing bar between the two white lines), lower the two sidebars that help to keep the bag correctly positioned.


In the picture below you can see an example of a bag correctly positioned.

Now the operator can insert the injection nozzle into the bag by rotating 90° clockwise until it stops and locks into position.

Check again that the open side of the bag is positioned on the sealing bar and make sure that the edge of the bag protrudes approximately 1 cm, as shown below.

Position the upper & lower sides of the bag over the sealing bar and insert your hand inside the bag to flatten the material under the nozzle against the metallic slope. It is very important that the 2 sides of the bags are not attached to each other and that the lower side is flat.

A.2. How to correctly seal a vacuum bag

The internal part of the bag should be clean and dry to achieve a good sealing. It is also very important to correctly position the bag over the sealing bar, as explained in this section.

At the end of a program, it is required to verify that the sealing is uniform all along the side of the bag.

In case the sealing of the bag is not correct:

- Open the bag by using some scissors.
- Clean the upper internal part of the bag as explained above.
- Run immediately a program "Vacuum Only" to seal again the bag (leaving the formalin inside the bag). Choose the default program P3. See the chapter 3.6 to run a "Vacuum Only" program.

B. APPENDIX B - "VACUUM ONLY" PROGRAMS SELECTION GUIDELINES**Human Organs:**

P1	P2	P3	P4
Kidney (P3/P4 if with cyst)	Placenta	Lung	Cysts > 10cm in diameter
Breast	Prostate	Brain	
Uterus	Intestine	Lymph node	
Liver (P3/P4 if with cyst)	Testicle	Closed colon	
Muscle	Appendix	Gallbladder	
Thyroid	Bone (not sharp)**	Eye	
Thymus	Stomach/ Esophagus	Cysts <10cm in diameter	
Skin	Ovary	Fetus	
Tongue	Tonsil	Scraping/TURP	
Spleen			
Open colon*			
Pancreas			
Bladder			
Heart			

* If the colon is clamped, please use the program P3.

** Remember to cover the sharpened edges (if present) with gauze to protect the vacuum bag.


If the nature of the specimen is uncertain, use the program P3.

This guide is intended for specimens directly packaged with Milestone vacuum sealed plastic bags.

See the chapter 3.6 to run a "Vacuum Only" program.

C. APPENDIX C – SPECIMENS TRANSPORTATION GUIDELINES

If the time required by the courier for the specimens' transportation to the Anatomic Pathology is:

- **less than 30 minutes:** put the specimens in the Thermal bag with the cool packs at +4°C ,
- **less than 1 hour:** put the specimens into the Refrigerator, close to the SealSAFE, for at least 2 Hr (pre-cooling procedure) and later place them in the Thermal bag with the cool packs at +4°C,
- **more than 3 hours:** use of the Mobile Refrigerator connected to an external power supply (12Vdc, 115/230V±10% ~).

With the system at controlled temperature designed by Milestone s.r.l. is ensured a perfect preservation of the morphology, antigenicity, DNA/RNA and proteins are preserved up to 48/72 hours in an optimal manner (depending on the initial conditions and the type of sample).

In none of the three cases described above the total time of the transport must exceed 72 hours (including the time spent in the Surgery Room/Refrigerator of the Laboratory of Histology). For proper preservation of the morphology, the samples must be packed and stored at 4°C ±2°C.

D. APPENDIX D – PRINTER AND BARCODE READER


D.1. How to use the printer

A built-in label printer is included for complete documentation of run by a self adhesive label that identifies the sample bag.

1D barcodes are printed by this printer.

A printer kit (CODE 68331/A) is included (10 label rolls); see the chapter 2.2 for order more.

The printer has an ABS casing with a front cover which opens to allow access to the paper roll and print head. The control panel is located on the front (see the figure below) and has a FEED key (3) and an OPEN key (1) for paper roll compartment opening. The OPEN key is made of transparent plastic material so it works also a STATUS Led (2).

- FEED button (3): When the LINE FEED is pressed, the printer advances the paper.
- OPEN button (1): press the OPEN key to open the paper roll compartment.
- STATUS led (2): displays printer hardware status. In case of malfunction. The flash frequency changes as shown in table below.

STATUS LED	COLOR	DESCRIPTION	
	OFF	Printer OFF	
	GREEN	Printer ON: no error	
	GREEN	COMMUNICATION STATUS	
		Flushing	Description
		Slow	Paper end
		Fast	Heating over temperature or incorrect power supply voltage

At the installation time, a roll is already located into the printer and 9 additional rolls are included separately in the shipment.

When a roll is over, it's required to replace it if you want to continue printing. See the following instructions on how to change the paper roll.

If you don't want to use the print function or if you finish the paper roll and you don't have another one to insert and so the printer will stay without any roll for a period of time, it's required to disable the printing function in Owner Settings (refer to the chapter 3.4 for instructions).

D.1.1. Storage of paper rolls

Roll of paper for the printer (component)

Transportation conditions

Joint transportation with the instrument is allowed

Temperature: -20°C to +60°C

Humidity: up to 50% (at 40°C) non-condensing

Storage conditions

Store separately from the instrument


Temperature: +20 to 40°C.

Humidity: up to 50% (at 40°C) non-condensing.

D.1.2. Changing the paper roll

To change paper roll, proceed as follows:

- Open the paper roll compartment by pressing the OPEN key.


- Move down the plastic cover to the maximum opening.


Before inserting a new roll of paper into the printer, cut the first part (20cm) that may be sticky.

- Place the paper roll making sure that it unrolls in the proper direction.

- Take out the paper and close the plastic cover.
- Push on the plastic cover to lock it.
- Tear off the exceeding paper using the jagged edge.

D.1.3. Printer troubleshooting

It's required to use the jagged edge when tearing off the paper.

In case the printer doesn't advance the paper correctly, verify that the roll has been inserted correctly by following the operations described in the previous chapter.

It is also suggested to press the button FEED of the printer after a sealing. The printer advances a little the paper and so at the moment of the following printing process, the printer should advance the paper correctly because it has been already advanced.

E. APPENDIX E – BAG SEALER WITH HANDLE

How to use the bag sealer with handle

The bag sealer with handle (CODE 68559) will be included in SealSAFE unit in quantity 2. If necessary, it will be possible to order it as accessory additionally, minimum order 10.


The bag sealer can be used to keep closed temporary a bag containing or not a specimen, dirty or not by fixative.

How to use the Bag sealer:


- Bend the upper part of the bag as shown into the picture.

- Insert the yellow spindle of the bag sealer into the bended part of the bag.

- Then push the bag sealer along the bended part of the bag.

- The bag sealer will keep closed all the opened side of the bag.


Below the final result.

It is possible to use the handle of the bag sealer to move a bag inside the laboratory but it is required to keep one hand below the base of the bag for safety.

F. APPENDIX F - RAPID CONNECTOR WITH HOSE SUITABLE FOR POLITAINER TANKS (CODE 68630)**Using of Rapid connector with hose suitable for Politainer tanks (CODE 68630)**

The rapid connector (CODE 68630) provided is a new fast coupling, to be used in the place of the default ones. See the image below.

By using this new accessory, the connection of formalin tanks with SealSAFE becomes safer because the operator has only to connect the tank to the rapid connector, without opening it. Therefore, formalin fumes are strongly reduced in the air during the tanks replacement procedure.

Formalin tanks are not provided.

The quick tank coupling is Cornelius Keg (with ball lock).

HOW TO CONNECT A TANK INTO THE DRAWER WHEN RAPID CONNECTOR IS USED:


- Connect the rapid connector in position 1 or 2 into the drawer, as shown below.

- Position the formalin tank keeping on the top the ball lock.
- Connect the rapid connector with the tank.

- Rotate the tank and position it as shown in the picture below otherwise the loading is not correct.

- When the tank is empty, rotate it in vertical position, unplug the rapid connector and position a new one following the steps above.

G. APPENDIX G – THE LIST OF NATIONAL STANDARDS APPLIED BY THE MANUFACTURER

The following standards apply:

- 98/79/EC In vitro diagnostic medical devices Directive
- 2014/35/EC Low Voltage Directive
- 2014/30/EC Electromagnetic compatibility (EMC) Directive
- 2011/65/EC Directive on the restriction of the use of certain hazardous substances in electrical and electronic equipment
- CEI EN 61010-1:2001 Safety requirements for electrical equipment for measurement, control and laboratory use. General requirements.
- CEI EN 61010-2-101:2002 Safety requirements for electrical equipment for measurement, control and laboratory use. Particular requirements for in vitro diagnostic (IVD) medical equipment.
- CEI EN 61326-1:2006 Electrical equipment for measurement, control and laboratory use. EMC requirements. Part 1. General requirements.
- CEI EN 61326-2-6:2006 Electrical equipment for measurement, control and laboratory use. EMC requirements. Part 2-6. Particular requirements. In vitro diagnostic (IVD) medical equipment.
- EN IEC 62304:2006 (/2013) Medical devices. Software. Life cycle process.
- IEC 60825-1-2013 Safety of laser products. Part 1. Equipment classification, requirements and user's guide
- GOST IEC 61010-1-2014 Safety requirements for electrical equipment for measurement, control, and laboratory use. Part 1. General requirements.
- GOST IEC 61010-2-101-2013 Safety requirements for electrical equipment for measurement, control, and laboratory use. Part 2-101. Particular requirements for in vitro diagnostic (IVD) medical equipment.
- GOST R IEC 61326-1-2014 Electrical equipment for measurement, control and laboratory use. EMC requirements. Part 1. General requirements.
- GOST R IEC 61326-2-6-2014 Electrical equipment for measurement, control and laboratory use. EMC requirements. Part 2-6. Particular requirements. In vitro diagnostic (IVD) medical equipment.
- GOST 14254 -2015 Degrees of protection provided by enclosures (IP Code).
- GOST 20790-93 Medical instruments, apparatus and equipment.

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