Siemens Healthcare GmbH
Henkestr. 127,
91052 Erlangen,
Germany

User Manual addendum for Russia

SOMATOM Drive Computed Tomography X-ray system with accessories.

Supplemental information regarding Instruction for use for Russia

1. The name of the Medical Device.

SOMATOM Drive Computed Tomography X-ray system with accessories.

 CAUTION	
Not observing the Instructions for Use of the system, system options and accessories!	
Injury to the patient.	
◆	Always use the Instructions for Use in conjunction with the Instructions for Use of the particular units used.
◆	Follow the safety instructions.

2. Manufacturing information of the Medical Device.

Legal manufacturer:

Siemens Healthcare GmbH
Henkestr. 127,
91052 Erlangen,
Germany
Web: <https://www.healthcare.siemens.com>

Manufacturing site:

Siemens Healthcare GmbH
Computed Tomography (CT)
Siemensstr. 1,
91301 Forchheim
Germany
Tel.: + 49 9191 18-0
Web: www.healthcare.siemens.com/computed-tomography

The name and address of the local representative in Russia:

Siemens Healthcare LLC
Dubinskaya street 96
115093 Moscow,
Russia
Tel.: +7 495 737 10 00
Fax: +7 495 737 13 20
<https://www.siemens-healthineers.com/ru/>

3. Designation of the Medical Device.

3.1. Intended use

This computed tomography system is intended to generate and process cross-sectional images of patients by computer reconstruction of x-ray transmission data.

The images delivered by the system can be used by a trained physician as an aid in diagnosis.

3.2. Indications for use

This computed tomography system is intended to generate and process cross-sectional images of patients by computer reconstruction of x-ray transmission data.

The images delivered by the system can be used by a trained physician as an aid in diagnosis.

The images delivered by the system can be used by trained staff as an aid in diagnosis, treatment preparation and radiation therapy planning.

This CT system can be used for low dose lung cancer screening in high risk populations.

3.3. Contraindications:

There are no known specific situations that contraindicate the use of the CT system.

3.4. Possible side effects:

There are no known side effects.

3.5. Conditions for use:

SOMATOM Drive must only be used by persons with the necessary specialist skills who have had appropriate application training, according to country-specific regulations, for example, physicians, trained radiologists, or trained technologists.

As the user, you must have the necessary qualifications (for example, training in a relevant medical profession). You must also have been instructed in the use of the SOMATOM system (for example, detailed knowledge of the instructions for use).

SOMATOM Drive can be operated in the X-ray department of a hospital or a doctor's surgery.

3.6. Patient population

No patient population restrictions are known.

3.7. Proper and improper use

Improper use means that the systems or parts of the systems are used against their intended use and indications for use.

 CAUTION
Impermissible or faulty manipulations or changes of the hardware or software, or connection of the system to a network can cause the system to malfunction!
Unauthorized access, injury to the patient or damage to the equipment.
◆ Make sure that all necessary precautions with respect to the existing level of security are considered when adding a functionality or altering the shipped configuration. ◆ Do not open or remove the cover of the equipment nor install third-party software.

Operation in patient vicinity

Only medical components of the system (gantry, patient table, power distribution cabinets (PDCs), accessories) are suitable for use within the patient environment.

Operation mode

The system is developed for continuous operation with intermittent loading. Thus, it should be switched off and restarted regularly. For details, see System Owner Manual.

3.8. Medical practice

The medical information presented in these instructions for use is for illustration only and is not intended to be used by the reader for the practice of medicine.

Any health care practitioner reading this medical information is reminded that they must use their own learning, training, and expertise in dealing with their individual patients. This material does not substitute for that duty and is not intended by Siemens to be used for any purpose in that regard.

The doses of the drugs mentioned are consistent with their approval labeling. The treating physician bears the sole responsibility for the diagnosis and treatment of patients, including drugs and doses prescribed in connection with such use.

4. Classification of the Medical Device.

4.1. Potential risk class:

Classification according to Annex IX, MDD

Class:	IIB
--------	-----

4.2. Class of the potential risk of the medical device application:

According to GOST R 31508-2012 Medical devices:

Device is intended for transmission of energy in the form of ionizing radiation, it belong to the class 2b.

4.3. Type of the Medical Device in accordance with the Russian Nomenclature classification:

135190

4.4. Class of the Medical Device in accordance with IEC 60601-1:

In according to IEC 60601-1: CLASS I.

4.5. Type of contact with the human body:

The Medical Device is non-sterile, nonexpendable.

Device surfaces that can come into contact with patients are evaluated and proven to be biocompatible according to ISO 10993.

4.6. Software classification

4.6.1 Software safety class

Class B according to IEC 62304

4.6.2 Software Version

VB20.

4.6.3 Date of release

05/2019.

4.6.4 Size in GB

50 GB

5. Operation (Principles of the MD operation)

This computed tomography system is intended to generate and process cross-sectional images of patients by computer reconstruction of x-ray transmission data.

The images delivered by the system can be used by a trained physician as an aid in diagnosis.

The images delivered by the system can be used by trained staff as an aid in diagnosis, treatment preparation and radiation therapy planning.


A CT examination consists of the following steps:

- Patient registration
- Selecting the scan protocol
- Checking the patient position
- Gantry operator panel
- Selecting the predefined patient table position
- Positioning the patient
- Loading the patient using the gantry control panel
- Using the 3D camera
- Setting the patient table position to zero
- Unloading the patient
- Checking the topogram parameters
- Loading the topogram
- Acquiring the topogram
- Planning a spiral scan
- Performing a spiral scan
- Reconstructing images
- Generating 3D reconstructions
- Ending the examination

Spiral scan

- Selecting the chronicle entry
- Moving the scan range
- Changing the length of the scan range or recon range
- Changing the field of view (FoV)
- Checking parameters
- Switching on API

6. The overview of the Medical Device.

- (1) Accessories;
- (2) Scan unit;
- (3) 3D camera;
- (4) Heat exchanger;
- (5) Power distribution cabinet;
- (6) Image reconstruction system;
- (7) Patient table;
- (8) Console.

Components or software functionalities that may not be part of your system configuration are nevertheless described in the present instructions for use. It may be that not all of them are marked explicitly as optional. The availability of these components or software functionalities depends on your purchase contract.

7. List of components and accessories of the Medical Device:

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			Height: 160 mm
25.	Опора для головы/предплечья.	Head-Arm-Support	Head-Arm-Support is used to place the hands above the patient's head in a cross position. Weight: 1,8 kg Width: 450 mm Length: 480 mm Height: 280 mm
26.	Комплект для фиксации рук пациента в положении за головой: - опора для головы/предплечья; - фиксирующие ремни, 2 шт.; - подушка.	Head-Arm-Support; Hand Fixation Straps for Head-Arm-Support; Cushion for Head-Arm-Support;	Are used to place the hands above the patient's head in a cross position. Head-Arm-Support: Weight: 1,5 kg Width: 480 mm Length: 455 mm Height: 255 mm Hand Fixation Straps for Head-Arm-Support: Weight: 0,2 kg Width: 75/50/24mm Length: 370 mm Cushion for Head-Arm-Support: Weight: 0,3 kg Width: 241mm Length: 210mm Height: 41/91mm
27.	Переходник опоры для головы/предплечья.	Adapter for Head-Arm Support	Is used for examinations with the head-arm support on table tops with a flat surface. Weight: 0,3 kg Width: 475mm Length: 230mm
28.	Дополнительные мониторы основной и/или сателлитной консоли (не более 10 шт.).	Monitor 19"/24" Monitor 19 DSC1913-D	Additional monitor for easier operation.
29.	Передвижная стойка для мониторов (не более 5 шт.).	Monitor Cart Monitor Cart Mono/Dual	Dual Monitor Cart is used for secure mounting and easy use of an additional monitor in the office for research. Monitor cart mono Weight: 26.0 kg Height: 1357 mm (with screen holder) Monitor cart dual Weight: 28.0 kg Height: 1357 mm (with screen holder) Monitor cart mono (depending on system configuration) Weight: max 20.0 kg Width: 480 mm Width with Monitor: 420 - 575 mm Length: 560 mm Height: 1350 - 1400 mm (with screen holder) Monitor cart dual (depending on system configuration)

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			FAST Direction – Verifies by using the camera image that the patient is positioned in the correct scan direction as compared to the selected protocols. If not in the appropriate orientation, the technician is provided a warning signal.
65.	Боковые рельсы стола пациента (пара).	Table Side Rails	Bilateral pleurical table side rails for convenient mounting of accessories.
66.	Кронштейн для сосуда инфузионной системы.	Infusion bottle holder	The infusion bottle holder is used for up to two infusion bottles, with a maximum load of 20 N, which is approximately 2 kg (4,4 lbs), per hook. Weight: 1.5 kg Height: 1010 mm Length: 250 mm Width: 110
67.	Ключ лицензионный для активации опции возвратно-поступательного движения стола для динамических исследований.	Adaptive 4D Spiral	Adaptive 4D Spiral - a unique 4D Spiral scan mode that enables the SOMATOM Drive to extend beyond restraints experienced when utilizing a static detector and allows for up to 48 cm (18.89") dynamic CT coverage. This enables use not only in perfusion but also for advanced 4D CT DSA evaluations.
68.	Ключ лицензионный для активации опции динамических исследований и перфузии миокарда.	Heart Perfusion Scanning Adaptive 4D Spiral	Adaptive 4D Spiral - a unique 4D Spiral scan mode that enables the SOMATOM Drive to extend beyond restraints experienced when utilizing a static detector and allows for up to 48 cm (18.89") dynamic CT coverage. This enables use not only in perfusion but also for advanced 4D CT DSA evaluations. Heart Perfusion Scanning - is a mode, for quantitative, dynamic myocardial stress perfusion imaging, is a sequence shuttle mode to dynamically cover up to approximately twice the detector width for myocardial perfusion studies with sufficient temporal resolution of 75 ms even for high heart rates. For a heart rate of 63 beats per minute or less every single heartbeat and for a heart rate of greater than 63 beats per minute every second heartbeat, images were acquired. This it provides sufficient temporal resolution even for high heart rates.
69.	Ключ лицензионный для активации пакета спектральной визуализации (не более 10 шт.).	syngo Dual Energy Scan with Tin Filter FAST DE Results FAST DE (DE WS4D)	Dual Energy scanning options including Tin Filter modes and the applications to introduce Dual Energy as part of your routine daily workflow. Includes FAST DE results and FAST DE with Workstream 4D. syngo Dual Energy Scan with Tin Filter - option allows the use of both SOMATOM Drive X-ray sources simultaneously at

			different energies, while the Tin Filter reduces dose and at the same time increases energy separation by blocking unnecessary parts of the energy spectrum. syngo Dual Energy offers the possibility to acquire two spiral data sets simultaneously from a single scan running the tubes at 80/140 kV or 100/140 kV. The results are two data sets with diverse information. FAST DE Results - enables a straight forward Dual Energy workflow. You can select dedicated Dual Energy applications at the AWP and they will be sent directly to the PACS without any interaction needed. FAST DE Results is as easy as selecting a recon job and will enhance you daily workflow significantly. FAST DE (DE WS4D) - is a 4D workflow for the Dual Energy data with direct generation of axial, sagittal, coronal, or double-oblique images from standard Dual Energy scanning protocols. The Advantage: the reduction or elimination of time consuming, error prone, manual reconstruction steps and a reduction of data volume up to a factor of 10, since virtually all diagnostic information is captured in 3D slices.
70.	Ключ лицензионный для активации пакета расширенной комплексной диагностики в кардиологии.	HeartView Cardio BestPhase Plus syngo Calcium Scoring CT FAST Phase	Cardiac scanning options to enable the simple to use, routine cardiac CTA and calcium scoring workflows, including beta blocker independent scanning, one heart beat scanning and flex scanning to enable functional imaging at low doses. Includes: Heart View, Cardio Best Phase Plus and syngo Calcium Scoring CT. Heart View - The option supports ultra-fast FlashSpiral scanning for maximum dose saving and scan times down to a quarter heart beat. ECG-synchronized Flash Cardio Sequence for dose efficient but versatile low dose cardiac imaging, including high heart rates and functional evaluation. CT images of the heart are acquired in defined phases of the cardiac cycle by prospective ECG-triggered sequence scanning and retrospective ECG-gated spiral scanning. The ECG signal used for gating the CT images is acquired by an integrated ECG device. The ECG signal is displayed on the gantry front cover and the scan interface. Dual Source acquisition mode with single segment reconstruction enables heart-rate

independent temporal resolution of 75 ms (factor 2 higher than single source acquisition with same parameter). Down to 32.5 ms temporal resolution combining HeartView Flash acquisition with robust 2-segment reconstruction.

With prospective ECG-triggered scanning, quick scans are triggered by ECG signals. The revolutionary FlashSpiral Cardio mode collects data projections of the entire heart in ultrafast 250 ms within a single diastolic phase. This performance is a direct result of having 2 X-ray tubes, simultaneously collecting information, combined with unprecedented table feeds above 400 mm/s. All the while a true temporal resolution of 75 ms is applied for each individual image.

The Flash Cardio Sequence mode introduces the Siemens-only dual-step pulsing, that maintains a low dose level during the systolic phase to calculate ejection fraction in addition to coronary imaging, therefore combining low dose coronary imaging with functional information,

Retrospective gating is based on a continuous spiral scan with simultaneous ECG recording. The cardio spiral reconstruction allows volume imaging in selectable phases of the cardiac cycle. With retrospective ECG-gated spiral scans the ECG signal can be edited for improved image quality in the case of severe arrhythmia. A dedicated "Preview" tool enables the planning of the volume reconstruction during an optimal cardiac phase on the basis of axial single slices.

Adaptive ECG-synchronized dose modulation (pulsing) allowing for optimal dose savings of up to 50% compared to single source CT scanners, with advanced irregular and ectopic heartbeat detection algorithm. The 4% MinDose algorithm lets the user save even more dose on the patient examination for dedicated coronary analysis. A special algorithm allows to decrease dose during ECG-Pulsing to 4% of the full tube output, thus decreasing dose about an additional 25%. Data evaluation is only possible in the full dose frame.

Cardiac Best Phase Plus - is a dedicated software which automatically calculates and detects the optimal phase for motion-less coronary visualization. The phase is defined as either end-systole, end-diastole or both

			timepoints and automatically reconstructed. Also includes the DirectViewing™ tool, used for real time navigation through full volumes of up to 24 heart phases by using an integrated fast 3D volume viewer, available both on the Examination and Recon subtask card. Furthermore it provides VRT visualization of the coronaries with heart isolation in up to 20 phases within 15 seconds. DirectViewing™ completes the workflow of Cardio BestPhase™ by giving you the flexibility to individually visualize phases for all coronary arteries. Ca Scoring #AWP scanning - for quick calcium scoring based risk assessment and coronary age calculation (for single user at Acquisition Workplace) FAST Phase - is an intuitive solution fully integrated in the cardiac workflow for helping to quickly set up complex and usually time-consuming cardiac procedures reliably and repeatedly. By measuring of heart rate and rhythm the system automatically chooses the most appropriate phase of the heart cycle to scan and later to reconstruct based on fully configurable and customizable look up tables. This reduces complexity and helps to standardize cardiac acquisitions by making them more robust and less user-dependent.
71.	Специализированный оловянный фильтр.	Tin Filter	Tin Filter reduces dose and at the same time increases energy separation by blocking unnecessary parts of the energy spectrum
72.	Ключ лицензионный активации пакета быстрого автоматизированного планирования топограммы, исследования и реконструкций.	FAST Topo FAST Planning FAST 3D Align	leading applications to make positioning simple for our customers. FAST Topo - enables faster scan speeds in topograms, which prevents breath-hold artifacts. It also has the potential to decrease the topogram dose. FAST Planning - assists the scan and reconstruction planning, based on a topogram, to provide an easier, faster and standardized workflow in CT scanning. FAST Planning features the selection of the anatomical region of interest from a list prospectively defined scan and reconstruction ranges, automatic detection of the scan region(s) of interest and proposal of corresponding scan range(s) in the topogram (in a narrow or wide lateral FoV), optimized FoV and automatic iso-center adaptation for Head scans.

			FAST 3D Align - automatically corrects misalignment of anatomic structures, organs of the patient. It aligns those to fit it to the selected reconstruction plane for a highly automated reconstruction workflow. Additionally it minimizes the black area in the image through automatically adjusts recon field of view selection. FAST 3D Align works in combination with Workstream 4D
73.	Ключ лицензионный для активации пакета комплексной визуализации и диагностики.	SureView Flash Spiral 2x Adaptive Dose Shield CARE DOSE 4D CARE kV CARE Child CARE Dashboard CARE Bolus DoseMAP FAST Adjust XXL Mode 2cm	Including: SureView - Siemens unique pitch independent dose solution Flash Spiral - Ultra-fast spiral scanning in Dual Source mode with up to 450 mm/s, allows for additional dose saving especially in ECG-triggered scans Dual Adaptive Dose Shields - world's first dynamic tube collimation that protects the patient from clinically irrelevant radiation in every spiral scan. Both tubes are equipped with an Adaptive Dose Shield allowing reduction of direct peripheral exposure in Spiral CT. CARE Dose 4D - uses at first an automated adjustment of the dose level depending on patient size based on the attenuation values obtained from the standard (singular) topogram along the patient z axis. In addition CARE Dose4D uses a real-time adaptation of the tube current during the scan based on the actual attenuation of the X-ray beam measured around the patient. It delivers significant x-ray dose reduction (up to 68 %) possible for all body regions scanned compared with standard sequence or spiral scanning; CARE kV - First automated, organ-sensitive voltage setting to improve image quality and contrast-to-noise-ratio while optimizing dose and potentially reducing it by up to 60%. CARE Child - Dedicated pediatric CT imaging, including 70 kV scan modes and specific CARE Dose4D curves and Pediatric Protocols - Special examination protocols with 80 kV and a large range of adjustable mAs values for optimum adaptation of the radiation exposure to the age and weight of the child to be examined. CARE Profile - Visualization of the dose distribution along the topogram prior to the

			scan CARE Dashboard - Visualization of activated dose reduction features and technologies for each scan range of an examination to analyze and manage the dose to be applied in the scan CARE Bolus - Operating mode for CM-enhancement triggered data acquisition. The objective is optimum utilization of the contrast medium bolus in its "plateau" phase in the target organ. This option has been especially adapted to the increased speed and timing requirements resulting from the multirow capability and faster rotation. The CM enhancement is observed via monitoring scans in a user-defined ROI with a trigger threshold. As soon as the enhancement reaches its predefined threshold, the spiral scan is triggered as quickly as possible. License for software use on one modality. Dose MAP - Siemens unique solution to support dose monitoring, protocol management and benchmarking FAST Adjust - assists the user to handle system settings in a fast and easy way by automatically solving of conflicts within user defined limits by one single click on the FAST Adjust button. The limits for scan time and tube current per scan are defined via the Scan Protocol Assistant. FAST Adjust offers an undo functionality to return to previously set values. XXL Mode 2cm - Utilize the power of both tube to allow for low kV scanning across a wider range of patient body sizes. ADMIRE - (Advanced Modeled Iterative REconstruction) offers significant dose reduction and image quality improvement, as well as an everyday suitability.
74.	Ключ лицензионный для активации расширенного пакета комплексной визуализации и диагностики.	DistinCT - Sigma High Power incl. High Power 70 High Power 80 Dual Power 4cm 10kV Steps Tin Filter X-CARE CARE Contrast III	Includes the DistinCT - Sigma High Power containing, High Power 70 and High Power 80, Dual Power 4cm and 10kV Steps, as well as Tin Filter scanning and CARE Contrast DistinCT - Sigma High Power - combines unique tube related features of the Straton MX Sigma. High Power 70 - 650 mA is produced at 70 kV High Power 80 - 750 mA is produced at 80 kV Scanning with low kV has the potential to lower the dose. Scanning with low kV makes it possible to reduce contrast media. Using less contrast

			media has a possitive impact on the patients health and reduces expences for the hospital. Low kV scanning has a beneficial impact on image impression: Scanning with 80 kV might increase the lesion conspicuity or facilitate radiation dose reduction or support IQ improvement at contrast-enhanced hepatic CT Dual Power 4cm - utilizing both tubes at the same time the system doubles its power whilst maintaining the speed of the system. These CTA protocols enable the use of low kV imaging in a wider variety of patient sizes. Post interventional scanning is often required for vascular patients whose bodies are under stress. Low kV scanning not only reduces the dose for the patient but can support to reduce contrast media for patients with post operational renal impairment, where CT is still indicated. 10kV steps - The kVs can be automatically pre-selected through CARE kV in finest 10 kV steps based on patient body habitus and examination type for lowest possible dose at constant signal to noise ratio (image quality). 10 kV steps allow a more distinct set of parameters are found for each patient. X-CARE - Partial scanning to reduce direct X-ray exposure for the most dose-sensitive body regions, e.g. the breasts, thyroid gland or eye-lens.
75.	Ключ лицензионный для активации пакета выполнения различных функциональных исследований.	syngo Dual Energy Scan with Tin Filter FAST DE Results FAST DE (DE WS4D) HeartView Cardio BestPhase Plus syngo Calcium Scoring CT FAST Phase Heart Perfusion Scanning Adaptive 4D Spiral	Innovating for the Therapy, as every patient is different and they have their own distinct treatment requirements and responses. With SOMATOM Drive we combine our market leading applications, Dual Energy, Cardiac and Dynamic Imaging, to make quantative scanning and analysis more accessable and simpler to use. Including : syngo Dual Energy Scan with Tin Filter - option allows the use of both SOMATOM Drive X-ray sources simultaneously at different energies, while the Tin Filter reduces dose and at the same time increases energy separation by blocking unnecessary parts of the energy spectrum. syngo Dual Energy offers the possibility to acquire two spiral data sets simultaneously from a single scan running the tubes at 80/140 kV or 100/140 kV. The results are two data sets with diverse information.

		FAST DE Results - enables a straight forward Dual Energy workflow. You can select dedicated Dual Energy applications at the AWP and they will be sent directly to the PACS without any interaction needed. FAST DE Results is as easy as selecting a recon job and will enhance you daily workflow significantly. FAST DE (DE WS4D) - is a 4D workflow for the Dual Energy data with direct generation of axial, sagittal, coronal, or double-oblique images from standard Dual Energy scanning protocols. The Advantage: the reduction or elimination of time consuming, error prone, manual reconstruction steps and a reduction of data volume up to a factor of 10, since virtually all diagnostic information is captured in 3D slices. Heart View - CT images of the heart are acquired in defined times of the cardiac cycle by prospective ECG-triggered sequence scanning and retrospective ECG-gated spiral scanning. The ECG signal used for gating the CT images is acquired by an integrated ECG device. The ECG signal is displayed on the gantry Touch Panels and the scan interface. Dual Source acquisition mode with single segment reconstruction enables heart-rate independent temporal resolution of 75 ms (factor 2 higher than single source acquisition with same parameter). Down to 32.5 ms temporal resolution combining HeartView acquisition with robust 2-segment reconstruction. The Adaptive Cardio Sequence mode introduces the Siemens-only dual-step pulsing, that maintains a low dose level during the systolic phase to calculate ejection fraction in addition to coronary imaging, therefore combining low dose coronary imaging with functional information. Retrospective gating is based on a continuous spiral scan with simultaneous ECG recording. The cardio spiral reconstruction allows volume imaging in selectable phases of the cardiac cycle. With retrospective ECG-gated spiral scans the ECG signal can be edited for improved image quality in the case of severe arrhythmia. A dedicated "Preview" tool enables the planning of the volume reconstruction during an optimal cardiac phase on the basis of axial single slices.
--	--	--

		Adaptive ECG-synchronized dose modulation (pulsing) allowing for optimal dose savings of up to 50% compared to single source CT scanners, with advanced irregular and ectopic heartbeat detection algorithm. The 4% MinDose algorithm lets the user save even more dose on the patient examination for dedicated coronary analysis. A special algorithm allows to decrease dose during ECG-Pulsing to 4% of the full tube output, thus decreasing dose about an additional 25%. Data evaluation is only possible in the full dose frame. Cardiac Best Phase Plus - is a dedicated software which automatically calculates and detects the optimal phase for motion-less coronary visualization. The phase is defined as either end-systole, end-diastole or both timepoints and automatically reconstructed. Also includes the DirectViewing™ tool, used for real time navigation through full volumes of up to 24 heart phases by using an integrated fast 3D volume viewer, available both on the Examination and Recon subtask card. Furthermore it provides VRT visualization of the coronaries with heart isolation in up to 20 phases within 15 seconds. DirectViewing™ completes the workflow of Cardio BestPhase™ by giving you the flexibility to individually visualize phases for all coronary arteries. Ca Scoring #AWP scanning - for quick calcium scoring based risk assessment and coronary age calculation (for single user at Acquisition Workplace) FAST Phase - is an intuitive solution fully integrated in the cardiac workflow for helping to quickly set up complex and usually time-consuming cardiac procedures reliably and repeatedly. By measuring of heart rate and rhythm the system automatically choses the most appropriate phase of the heart cycle to scan and later to reconstruct based on fully configurable and customizable look up tables. This reduces complexity and helps to standardize cardiac acquisitions by making them more robust and less user-dependent. Adaptive 4D Spiral - a unique 4D Spiral scan mode that enables the SOMATOM Drive to extend beyond restraints experienced when utilizing a static detector and allows for up to 48 cm (18.89") dynamic CT coverage. This enables use not only in perfusion but also for
--	--	---

			advanced 4D CT DSA evaluations. Heart Perfusion Scanning - is a mode, for quantitative, dynamic myocardial stress perfusion imaging, is a sequence shuttle mode to dynamically cover up to approximately twice the detector width for myocardial perfusion studies with sufficient temporal resolution of 75 ms even for high heart rates. For a heart rate of 63 beats per minute or less every single heartbeat and for a heart rate of greater than 63 beats per minute every second heartbeat, images were acquired. This it provides sufficient temporal resolution even for high heart rates.
76.	Ключ лицензионный для активации пакета расширения режимов работы рентгеновских трубок.	High Power 70 High Power 80 Dual Power 4cm 10kV Steps	Sigma High Power enables the two Straton MX Sigma tubes to provide the full scope of the tube power. Includes: 10kV steps - The kVs can be automatically pre-selected through CARE kV in finest 10 kV steps based on patient body habitus and examination type for lowest possible dose at constant signal to noise ratio (image quality). 10 kV steps allow a more distinct set of parameters are found for each patient. Single energy voltage options of 70, 80, 90, 100, 110, 120, 130 and 140kV even in High Pitch Modes. High Power 70 - At 70kV the tube produces 650mA which is an increase on the original Straton tube of 150mA. This additional power will enable the use of low kVs in a greater number of patients, further increasing the utilization of dose and contrast media effective low kV scanning. High Power 80 - At 80kV the tube produces 750mA which is an increase on the original Straton tube of 150mA. This additional power will enable the use of low kVs in a greater number of patients, further increasing the utilization of dose and contrast media effective low kV scanning. Scanning with low kV has the potential to lower the dose. Scanning with low kV supports in the reduce contrast media. Using less contrast media may have a possitive impact on the patient's health and may provide a reduced expence for the hospital. Low kV scanning has a beneficial impact on image impression: Scanning with 80 kV might increase the lesion

			conspicuity or facilitate radiation dose reduction or support IQ improvement at contrast-enhanced hepatic CT Dual Power 4cm - utilizing both tubes at the same time the system doubles its power whilst maintaining the speed of the system. These CTA protocols enable the use of low kV imaging in a wider variety of patient sizes. Post interventional scanning is often required for vascular patients whose bodies are under stress. Low kV scanning not only reduces the dose for the patient but can support to reduce contrast media for patients with post operational renal impairment, where CT is still indicated. Sigma High Power also enables Tin Filter sanning.
77.	Ключ лицензионный для активации пакета расширения режимов работы генераторов.	120kV Tuning Dual Power 2cm 20kV Steps	With Sigma Power, the Somatom Drive comes equipped with the two Straton MX Sigma tubes, with the full scope of the tube power restricted. Includes: 20kV steps - The kVs can be automatically pre-selected through CARE kV with 20 kV steps based on patient body habitus and examination type for lowest possible dose at constant signal to noise ratio (image quality). Single energy voltage options of 70, 80, 100, 120 and 140kV even in High Pitch Modes. 120kV Tuning - The system has the maximum mA set at 120 kV - 833mA. 70kV is at 500mA 80kV is at 650mA
78.	Ключ лицензионный для активации пакета дополнительных инструментов для работы с изображениями.	VRT WorkStream 4D Extended Field of View	Innovating for the Radiologist, as every patient is different, Radiologists also vary, and they have their own distinct reconstruction and imaging requirements. With the SOMATOM Drive we combine our market leading applications to make reporting consistant, fast and simple for our customers Includes: VRT - Volume Rendering Technique standard on the 3D card. Workstream 4D - the standardized workflow guide for confident patient management. Up to 8 pre definable axial, coronal, sagittal and oblique MPR and MIP up to sub mm recon jobs possible. The Asynchronous Recon allows for multiple image reconstructions and reformats, parallel to scanning. With this feature, up to eight reconstruction job requests can be loaded into a scan protocol.

			Immediately upon completion of the scan acquisition, these reconstruction jobs are automatically executed in the background without delaying the start of next patient examination. WorkStream4D eliminates manual reconstruction steps and reduces the data volume up to a factor of 10, since virtually all diagnostic information is captured in 3D slices. Extended FoV - allows you to capture more information in just one exam, saving valuable time with emergency patients. Up to 78 cm field of view.
79.	Ключ лицензионный для активации расширенного пакета дополнительных инструментов для работы с изображениями.	FAST Spine #AWP iMAR #AWP HD FoV Pro #AWP z-UHR incl. UHR	advanced applications to make reporting of complex and atypical anatomical structures faster and simpler. Includes: iMAR #AWP The iMAR metal artifact reduction algorithm combines three successful approaches (beam hardening correction, normalized sinogram inpainting and frequency split). This allows to reduce metal artifacts caused by metal implants such as coils, metal screws and plates, dental fillings or implants. iMAR is compatible with extended FoV, the extended CT scale as well as the newest dose reduction feature. Along with the new algorithm comes the simple user interface of iMAR enabling easy reconstruction of clinical images with reduced metal artifacts. The iMAR metal artifact reduction algorithm combines three successful approaches to reduce metal artifacts: beam hardening correction (in sinogram regions of less severe metal attenuation), normalized sinogram inpainting (in sinogram regions of high metal attenuation), and frequency split (to mix back noise texture and sharp details that are potentially lost during inpainting). The correction process is then iteratively refined by repeating the normalized sinogram inpainting and the mixing steps thanks to the Adaptive Sinogram Mixing. Along with the new algorithm comes the simple user interface of iMAR. Besides the typical reconstruction parameters it only requires to select the desired protocol from a drop down menu which contains the following type of implants: dental fillings, neuro coil, thoracic coil, hip implants,

		extremity implants, pacemakers, spine implants and shoulder implants. FAST Spine #AWP Accurate and anatomically aligned preparation of spine recons with just a single click. FAST Spine provides various modes that automatically create anatomically orientated spine reconstructions based on a 3D volume. It provides an easier, faster and standardized workflow in CT scanning. FAST Spine features automatic segmentation of the spinal canal, automatic labeling of the vertebrae, anatomically oriented slices - (orthogonal to the spinal canal), coronal and sagittal reconstructions which refer to the curvature of the spinal column and more. All modes offer the possibility to adapt the results manually. FAST Spine requires Workstream 4D. z-UHR incl. UHR Siemens proprietary z-UHR (z-Ultra High Resolution) option enables previously unachievable image detail with an isotropic resolution of 30 lp/cm (0.17mm) at 0% MTF (+/- 10%). In addition, z-UHR Catphan® measurements demonstrate the industry's highest visible high-contrast resolution of 0.24 mm x 0.24 mm x 0.24 mm. The combination of z-Sharp Technology and z-UHR offers, in daily clinical routine, an isotropic detail in the range of research CsI-aSi flat-panel and Micro CT technology. z-UHR is intended for ultra-high resolution bone-imaging, in particular for small structures such as inner ear, joints or fractures of the bone. The option includes the z-UHR software as well as a dedicated tantalum detector comb.
--	--	---

8. Additional technical requirements

Parameter	SOMATOM Drive
Number of slices	2 x 128
Max. number of reconstructed slices	2 x 384
Maximum scanning field (cm)	50

Maximum reconstruction field (cm)	78
Cooling rate, MHU/min	7.3
Contrast difference	3 HU
2% MTF	UHR: - x/y-plane: 2% MTF 24.0 пар лин./см ($\pm 10\%$); - z-direction: 2% MTF 24.3 пар лин./см ($\pm 10\%$). Standard: - x/y-plane: 2% MTF 16.4 пар лин./см ($\pm 10\%$); - z-direction: 2% MTF 16.9 пар лин./см ($\pm 10\%$).
Number of projections	1152, 2304, 4608

9. Capabilities of the Medical Device with other medical devices.


List of accessory products approved by Accessory Solutions for CT Systems

Only qualified accessories are allowed to be operated in combination with Siemens Healthcare modalities. If the accessories, listed in this document are offered in combination with Siemens modalities, an "Article 12 EC Declaration" is mandatory.


	Supplier
Accutron CT-D Power Supply Version Accutron CT-D pedestal Accutron CT-D Ceiling Accutron CT-D Interface	Medtron
Stellant with Certegra WS pedestal Stellant with Certegra WS ceiling Interface: ISlgoo Centargo, pedestal - battery	Bayer (Medrad)
Optivantage DH pedestal Optivantage DH ceiling	Guerbet (Mallinckrodt)
Empower CTA+ pedestal Empower CTA+ ceiling CT Exprès with Interface, pedestal CT Exprès with Interface, ceiling mounted	Bracco
CT motion pedestal CT motion ceiling Interface	Ullrich
Dual Shot Alpha 7 Dual Shot Alpha 7s Dual Shot Alpha B200 Dual Shot GX Interface	Nemoto
Hardware Respiratory Gating System, international, AZ-733VI	Anzai
RGSC Varian Respiratory Gating System	Varian
Type-S Extension PosiFix Extension	Civio

10. Labeling Information and requirements:

Front side:


	
	(1) Follow the Instruction for Use. The user must consult and follow the Instruction for Use before using the CT system. An electronic copy of the Instruction for Use is available at http://www.siemens.com/medImaging-Manuals

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SIEMENS XXXXXXXX  (01) xxxxxxxxxxxxxxxxGTIN (240) xxxxxxxxMODEL (422) xxxORIGIN YY (21) xxxxxSERIAL (20) xxREV IVK Siemens Healthcare GmbH, Henkestr. 127, 91052 Erlangen, Germany Made in Germany	(4) Product identification label
	(5) Nationally Recognized Testing Laboratory (NRTL)
CE  CMII ID: XXXXX  Z XXX  N XXX  R xxx WW xxxxxxxx Contains FCC ID: xxxxxxxx IC: xxxxxxxx Complies with IMDA Standards DA xxxxx	(6) Radio frequency information Only for systems with an optional wireless i-Control.

Gantry frame:
The labels at the gantry frame are located beneath the gantry cover.

	(1) Product status label Refer to the original label that is applied on the device for information about the current product version. This label is an example only.
--	--

Rev.	00	01	02	03	04	05	06	07	08	09	10	11
12	13	14	15	16	17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32	33	34	35	36	37

 XXXXXXXX (01) xxxxxxxxxxxxxxxx GTIN (240) xxxxxxxx MODEL (422) xxx ORIGIN YY (21) xxxxx SERIAL (20) xx REV Siemens Healthcare GmbH, Henkestr. 127, 91052 Erlangen, Germany Made in Germany	(2) Product identification label
	(3) Approval label (USA and Canada, NRTL) NRTL – Nationally Recognized Testing Laboratory
 SOMATOM XXXXXXXXX (01) xxxxxxxxxxxxxxxx GTIN (240) xxxxxxxx MODEL (422) xxx ORIGIN YY (21) xxxxx SERIAL (20) xx REV Siemens Healthcare GmbH, Henkestr. 127, 91052 Erlangen, Germany Made in Germany	(4) Product identification label (System)


Tube housing:

The diagram shows a side view of the STRATON XX X device. Six labels are indicated by numbered circles and leader lines:


- (1) Product status label: Located at the top of the device.
- (2) Focal spot size label: Located on the upper left side of the device.
- (3) Anode/filter label: Located on the upper right side of the device.
- (4) kVp label: Located on the lower right side of the device.
- (5) Exposure time label: Located on the lower left side of the device.
- (6) mAs label: Located on the front of the device, below the exposure time label.

	(3) Product identification label
This product complies with DHHS regulations 21 CFR Subchapter J, applicable at date of manufacture. Manufactured: XXXXXXXXXXXX Siemens Healthcare GmbH Henkestr. 127, 91052 Erlangen Germany	(4) Approval label (USA, DHHS) DHHS – Department of Health and Human Services
	(5) Approval label (USA and Canada, NRTL) NRTL – Nationally Recognized Testing Laboratory
	(6) Product identification label

Tube collimator:

																																								
	(1) RoHS label																																							
This product complies with DHHS regulations 21 CFR Subchapter J, applicable at date of manufacture. Manufactured: XXXXXXXXXXXX Siemens Healthcare GmbH Henkestr. 127, 91052 Erlangen Germany	(2) Approval label (USA, DHHS) DHHS – Department of Health and Human Services																																							
X-Ray Filters fix: XX: X,X mm XX: X,X mm XX: X,X mm switchable: XX: X,X mm	(3) Filtering label Refer to the original label that is applied on the device for information about the X-ray filters. This label is an example only.																																							
SIEMENS XXXXXXXX  (01) xxxxxxxxxxxxxxxxGTIN (240) xxxxxxxx (422) xxx (21) xxxxx (20) xxx MODEL ORIGIN YY SERIAL REV IVK Siemens Healthcare GmbH Henkestr. 127, 91052 Erlangen, Germany Made in Germany	(4) Product identification label																																							
<table><tr><td>Rev.</td><td>00</td><td>01</td><td>02</td><td>03</td><td>04</td><td>05</td><td>06</td><td>07</td><td>08</td><td>09</td><td>10</td><td>11</td></tr><tr><td>12</td><td>13</td><td>14</td><td>15</td><td>16</td><td>17</td><td>18</td><td>19</td><td>20</td><td>21</td><td>22</td><td>23</td><td>24</td></tr><tr><td>25</td><td>26</td><td>27</td><td>28</td><td>29</td><td>30</td><td>31</td><td>32</td><td>33</td><td>34</td><td>35</td><td>36</td><td>37</td></tr></table>	Rev.	00	01	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	(5) Product status label Refer to the original label that is applied on the device for information about the current product version. This label is an example only.
Rev.	00	01	02	03	04	05	06	07	08	09	10	11																												
12	13	14	15	16	17	18	19	20	21	22	23	24																												
25	26	27	28	29	30	31	32	33	34	35	36	37																												

High voltage generator:

Rev.	00	01	02	03	04	05	06	07	08	09	10	11
12	13	14	15	16	17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32	33	34	35	36	37

(1) Product status label

Refer to the original label that is applied on the device for information about current product version. This label is an example only

This product complies with DHHS regulations 21 CFR Subchapter J, applicable at date of manufacture.
Manufactured: XXXXXXXXXXXX
Siemens Healthcare GmbH
Henkestr. 127, 91052 Erlangen
Germany

(2) Approval label (USA, DHHS)

DHHS – Department of Health and Human Services



XXXXXXXXXXXX
 (01) xxxxxxxxxxxxxxxx GTIN
(240) xxxxxxxxxxxx MODEL
(422) xxx ORIGIN YY
(21) xxx SERIAL
IVK

(3) Product identification label

Service Part xxxxxxxxxxxx
 (01) xxxxxxxxxxxxxxxx GTIN
(240) xxxxxxxxxxxx MODEL
(422) xxx ORIGIN YY
(21) xxx SERIAL



(4) Spare part identification label

Detector DMS A:

(1) Product identification label


Detector DMS B:

(1) Product identification label



Patient table:


Patient table with 2000 mm scan range


(1) Warning, crushing hazard for hand


Increased risk that hands can get pinched between table top and patient table





	(2) Caution, leg injuries Contusion of legs is possible when lowering the patient table. Caution: Consult accompanying documents.
	(3) Approval label (USA and Canada, NRTL) NRTL – Nationally Recognized Testing Laboratory
	(4) Duty cycle label The cyclic factor (c.d.f) specifies the ratio of the maximum operation time of the component to the total operation time of the system
This product complies with DHHS regulations 21 CFR Subchapter J, applicable at date of manufacture. Manufactured: XXXXXXXXXXXX Siemens Healthcare GmbH Henkestr. 127, 91052 Erlangen Germany	(5) Approval label (USA, DHHS) DHHS – Department of Health and Human Services
SIEMENS XXXXXXXX (01) xxxxxxxxxxxxxxxx GTIN (240) xxxxxxxx MODEL (422) xxx ORIGIN YY (21) xxxxx SERIAL IVK Siemens Healthcare GmbH Henkestr. 127, 91052 Erlangen, Germany Made in Germany	(6) Product identification label
	(7) Type B applied part Equipment is not suitable for direct cardiac application. Indicates the degree patient protection of this applied part (Type B applied part according to IEC 60601-1).

	(8) Product identification label
	(9) Total working load label The maximum permissible working load must not exceed 227 kg (500 lbs), which is approximately 226 N.
	(10) Equipotentiality label Indicates the location of a pin that can be used for equalization with other medical devices in accordance with IEC 60601-1 and DIN 42801-2.


Multi purpose table

	(1) Warning, crushing hazard for hand Increased risk that hands can get pinched between table top and patient table
	(2) Product identification label
	(3) Equipotentiality label Indicates the location of a pin that can be used for equalization with other medical devices in accordance with IEC 60601-1 and DIN 42801-2.
-or-  -or- 	(4) Total working load label The maximum permissible working load must not exceed 227 kg (500 lbs), which is approximately 2226 N. -or- The maximum permissible working load must not exceed 307 kg (667 lbs), which is approximately 3010 N.
	(5) Type B applied part Equipment is not suitable for direct cardiac application. Indicates the degree of patient protection of this applied part (Type B applied part according to IEC 60601-1).

	(6) Approval label (USA and Canada, NRTL) NRTL – Nationally Recognized Testing Laboratory
	(7) Duty cycle label The cyclic duration factor (c.d.f) specifies the ratio of the maximum operation time of the components to the total operation time of the system.
	(8) Product identification label
	(9) Approval label (USA, DHHS) DHHS – Department of Health and Human Services
	(10) Caution, leg injuries Contusion of legs is possible when lowering the patient table. Caution: Consult accompanying documents.

Physiological measurement module:

(1) Warning label

This equipment may be dangerous to patient and user unless safe expose factors, Instructions for Use, and maintenance schedules are observed.

The Type BF applied part is not defibrillation proof.

(2) ECG socket

(3) Type BF applied part

Equipment is not suitable for direct cardiac application. Indicates the degree of patient protection of this applied part (Type BF applied part according to IEC 60601-1).


(4) ECG channel selection

(5) Equipotentially label

The label indicates the location of a pin for potential equalization with other medical devices in accordance with IEC 60601-1 and DIN 42801-2.

Physiological measurement module for multi-purpose patient table:

(1) Warning label

This equipment may be dangerous to patient and user unless safe expose factors, Instructions for Use, and maintenance schedules are observed.

The Type BF applied part is not defibrillation proof.


(2) ECG socket

(3) Type BF applied part

Equipment is not suitable for direct cardiac application. Indicates the degree of patient protection of this applied part (Type BF applied part according to IEC 60601-1).

	 (4) ECG channel selection
	(5) Equipotentially label  The label indicates the location of a pin for potential equalization with other medical devices in accordance with IEC 60601-1 and DIN 42801-2.
	(6) Warning, crushing hazard for hand  Increased risk that hands can get pinched between table top and patient table.

11.Artwork of the MD labeling in Russian

Additionally, a Russian label is applied with the following content (please note that Labeling mentioned below is applied only for the SOMATOM Drive systems which are delivered to Russia not earlier than April 2019):

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12. Environmental Protection Requirements:

The CT system and its components comply with the following product environmental protection regulations.

Statutory regulation	Title
2006/66/EC:2006	Directive 2006/66/EC of the European Parliament and of the Council of 6 September 2006 on batteries and accumulators and waste batteries and accumulators and repealing Directive 91/157/EEC
2011/65/EU:2011	Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment
1907/2006/EC:2006	Regulation (EC) No 1907/2006 of the European Parliament and of the Council of 18 December 2006 concerning the Registration, Evaluation, Authorization and Restriction of Chemicals (REACH), establishing a European Chemicals Agency, amending Directive 1999/45/EC and repealing Council Regulation (EEC) No 793/93 and Commission Regulation (EC) No 1488/94 as well as Council Directive 76/769/EEC and Commission Directives 91/155/EEC, 93/67/EEC, 93/105/EC and 2000/21/EC

For compliance with existing laws and requirements regarding product related environmental protection is a basic requirement. An objective of the Healthcare policy involves preventing or minimizing environmental pollution that might be caused by production, use and end-of life management of Healthcare products. Consistent adherence to measures of product related environmental protection may be key to overall cost savings over the life of Healthcare products.

The requirements described in this document are aimed at achieving an efficient implementation of product related environmental protection regulations. These requirements specifically apply to the following:

- Environmental requirements regarding products characteristics;
- The identification, restriction, and documentation of materials used in products;
- Specific requirements for the suppliers;
- Requirements for supplemental accompanying documents with respect to environmental information.

Product environmental performance:

Product related environmental improvement targets. The Organizational Unit shall inventory and rank product environmental aspects according to environmental significance to the extent that data is available. The ranking shall be performed for each standard configured model (i.e., without regard to further options). The Organizational Unit shall document the inventory and ranking of product environmental aspects.

The Organizational Unit shall select and document at least one formal environmental improvement target prior to the development phase, from the existing or similar model to the new model based on aspect significance, technical and economic feasibility.

The aspect inventory and ranking shall be checked and modified if there are significant changes in the environmental aspects when:

- Product specifications change hardware or
- An existing or new supplier changes the characteristics of a part (system component).

Environmental improvement targets and target implementation plans shall be entered into the Siemens Healthcare EHS management system.

Risk analysis:

Environmentally relevant risks, such as from emissions under normal use, radiation generation or leakage, or caused by the use of critical materials, or risks occurring during end of life treatment shall be considered in the scope of risk analysis. This risk analysis shall consider at least the production, use and end of life unit environmental life cycle phases.

Exclusion of restricted materials in prohibited applications:

Defined restricted materials that are prohibited in specified applications above specified thresholds shall not be used in products. If a restricted materials is found to be used in a product for a prohibited application, the restricted materials must be eliminated from this product before it may be sold into a market that prohibits its use.

To ensure market access, all newly designed parts and components where the release of functional specification is after 2008-04-01 shall adhere to the RoHS substance restrictions, except for existent exemptions or exemptions in future are in place.

Use of critical materials in non-prohibited applications:

Defined critical materials specified in Attachment A1 and in ProMed tool Group 9.1 should not be used in products if a non-critical materials is available that meets the technical specifications at equal or better performance and equal or lower total cost. If defined critical materials are used, the type, amount, and location of defined critical substances in the product and their application shall be recorded and maintained.

Recording of product and packaging ingredients:

The Organizational Unit shall record and document the type and amount of any defined critical substances and selective treatment parts used as accurately as possible. Also, it shall obtain information regarding all noncritical materials but making up at least 90 % by weight of delivered product and at least 90% by weight of its packaging (Separate lists for product and packaging ingredients is required). If this information cannot be obtained after reasonable efforts, the Organizational Unit shall estimate the missing information.

This information shall be documented using the ProMed tool at least for each system component making up a product, except for system components without critical materials, without selective treatment parts and with negligible mass compared to the product. To achieve this, it may be necessary to assess subparts of system components too. If a part is purchased or supplied by another company or facility, the material information should be for the part as supplied rather than for its subparts and components.

Product ingredients of system components shall be collated for the complete product in a typical configuration.

The information about product and packaging ingredients shall be maintained using the ProMed tool.

The Organizational Unit shall document the net mass for each system component making up a product and for non-local trade goods.

The Organizational Unit shall document type and mass of packaging for each system component making up a product and for non-local trade goods. Excluded is packaging that always comes back to the production facilities.

Use of disposables and consumables:

Disposables and consumables required for the operation of products must be minimized to the extent possible. If used, the type and estimated period of usage of disposables and consumables shall be identified and available to users and/or the customer service.

Use of batteries:

The use of restricted batteries has been avoided.

The product does not contain:

- batteries or accumulators containing more than 0.002% (20 ppm) of cadmium by weight except for batteries or accumulators that are intended to be used in products that are included in a Medical Device Directive(MDD) CE declaration, or
- batteries or accumulators containing more than 0.0005% (5 ppm) of mercury by weight except for button cells that shall not contain more than 2% (20000 ppm) of mercury by weight.

The Organizational Unit shall document type and weight of batteries and accumulators for each system component making up a product and for non-local trade goods.

Labeling of plastic parts:

Thermoplastic parts shall be labeled per ISO 11469 if economically feasible (especially when molded). Details of the labeling, such as the threshold mass of parts to be labeled, size of the labels, etc. shall be defined by the Organizational Unit.

Product labeling:

Product is labeled with the symbol "crossed-out wheeled bin" in combination with the solid bar as specified in EN 50419. The symbol may be printed on product packaging or included with user documentation if the product size or function does not practically allow for labeling on the product itself.

In addition the product shall be labeled with the identification of the producer. This can be done by using the product identification label (System IVK Label) or by using the "SIEMENS" brand name.

Environmental relevant experiences acquired in the product use:

The Service organization shall provide feedback to manufacturing Organizational Units on environmental performance serviceability issues using the installed product base tool or other data from the field. During the development of new products, environmentally relevant experiences based on use and end of life unit environmental life cycle phases of current generation products and comparable products should be taken into consideration by reviewing service metrics in the installed product base tool, complaints or other field data.

13. Sterility:

The medical device is non-sterile.

14. List of equipment and measurement tools necessary for quality control of the medical device:

Accessories for quality measurement:

You require the set of phantoms for the quality measurement and possibly additional accessories for the constancy test.

Accessories for daily quality measurement:

The following components are required for the daily quality measurement:

- Phantom holder
- Set of phantoms

Accessories for the monthly constancy test:

For the monthly constancy test, you require the following phantoms and accessories:

- Phantom holder;
- Set of phantoms;
- Patient equivalent load not exceeding 135 kg (297 lbs) (a weight of 100 kg (220 lbs) is recommended by Siemens);
- Ruler, 40 cm or longer.

Accessories for the Low Contrast test

For the low contrast test, you require the following phantoms and accessories:

- Phantom holder
- Low contrast phantom


15. Transportation conditions and requirements.

15.1. ACCELERATION

To avoid system damage, the following conditions are mandatory for transport:

- Maximum acceleration according to DIN EN 60721-3-2, Class 2M2.

Directions of impact (example):

(X) Permissible impact: 4 g (corresponds to shunting a train)

(Z) Permissible impact: 2 g

Transport all CT components in the respective transport packaging provided by the manufacturer. When transporting the CT system across country terrain, only use trucks with air shock suspension. Additionally, observe the "Safety information" section in the first chapter "General information".

15.2. CLIMATIC CONDITIONS FOR TRANSPORT AND STORAGE

Standard climatic conditions

The climatic requirements for storing the CT system are identical to the climatic requirements for transportation.

The data for transporting and storing the CT system are applicable only if the system is shipped free of damage in the transport packaging provided by the manufacturer.

Tab. 58 Climatic conditions: Transport and storage


Designation	Range
Temperature	-20 °C to +50 °C
Relative humidity	10% to 90%
Temperature gradient	10 K/h
Barometric pressure	700 hPa to 1060 hPa
Max. storage time	2 months

All components have to be stored in the respective transport packaging provided by the manufacturer.

Use equivalent packaging for intermediate storage if the original packaging provided by the manufacturer is no longer available.

Fig. 72: Transport and storage conditions (climatogram)

(1) 60 g/m³ absolute humidity

After moving the CT system from a cold environment into a warm room, make sure that the entire CT system is dry before switching it ON.

Extended climatic conditions:

If the temperature range during transport and storage exceeds the specified values:

- from -20 °C to -45 °C and
- from 50 °C to 70 °C

special packaging must be used. This limits the temperature for the system to the specified values. Verify that the components to be transported can withstand these temperatures without damage.

15.3. TRANSPORTING CT SYSTEM COMPONENTS

The project manager is responsible for an unobstructed move-in of the CT system components and the availability of the necessary transport equipment.

The project manager determines:

- If a crane is necessary for unloading the gantry from the truck;
- If the gantry transporting device is available;
- the delivery route into the examination room;
- If the floor protection is necessary (e.g., metal plates), and
- If doors, corridors, and elevators are suitable for the delivery of the CT system components.

Designation	Description
Gantry and patient table	• Crane suspension for unloading the gantry from a truck; • Gantry with a special transport device; • Patient table with premounted transport rollers at the table base.
Other CT system components	The dimensions, of the LCB, cabinets, operating console, components of the imaging system, tower PCs, monitors and keyboards do not present difficulties.

15.4. TRANSPORTING THE GANTRY

Tab. 60 Gantry transport data


Item	Description
Safety information	Observe the "Safety Information" section in the first chapter, "General Information".
Transport device	A special transport device is available and orderable for transporting the gantry. Refer to the CT Basic PG, print number CT00-000.891.04 for details regarding: ■ Order procurement ■ Functions

Item	Description
Safety information	Observe the "Safety Information" section in the first chapter, "General Information".
Transport dimensions	The following list shows the actual overall transport dimensions of the gantry (with covers) with the transport device lowered to the minimum ground clearance of 7 mm: - Standard transport (rollers are swiveled out) NOTE: The touch panels on the front and rear of the gantry do not affect the total transport width. 2645 x 1801 x 1999 mm (L x W x H) - Transport through narrow passages (rollers are swiveled in): NOTE: The touch panels on the front and rear side of the gantry affect the total transport width. a) With touch panels on the front side of the gantry only: 3331 x 1220 x 1999 mm (L x W x H) b) With touch panels on the front and rear side of the gantry: 3331 x 1256 x 1999 mm (L x W x H) NOTES: - Based on the aforementioned transport dimensions, the project manager must check all of the inner widths and heights in the building for moving in the gantry. - The transport device with the gantry can be lowered to a minimum ground clearance of 7 mm.
Transport weights	- Gantry with the transport device (without packaging): < 2865 kg - Transport device: < 285 kg
Floor load during transport	The maximum possible floor load (two-point load) per roller that may arise during gantry transport is 16 kN. NOTE: If necessary, cover the transport route with metal sheets to distribute the load.

Item	Description
Safety information	Observe the "Safety Information" section in the first chapter, "General Information".
Reducing the gantry weight	In special cases it is required to reduce the gantry weight by a few kilograms, e.g. to use a freight elevator with a limited freight load. This can be done by removing certain parts of the gantry covers. The following list shows the weight of the removable covers: ■ Hood, right: 15.76 kg ■ Hood, left: 12.71 kg ■ Side cover, left: 8.67 kg ■ Side cover, right: 10.95 kg ■ Front cover, right: 8.21 kg ■ Rear cover, right: 10.64 kg ■ Funnel: 13.36 kg ■ Baseboard cover set: 13.79 kg NOTE: A service engineer is required for disassembly and reassembly. Store the removed covers carefully at a safe location to avoid damage during transport of the gantry.

Item	Description
Safety information	Observe the "Safety Information" section in the first chapter, "General Information".

Fig. 73: Removing the gantry covers: Cover designations (gantry design with touch panels)

- (A) Front
(B) Back
(1) Hood right
(2) Hood left
(3) Side cover left
(4) Side cover right
(5) Front cover right
(6) Rear cover right
(7) Funnel
(8) Baseboard cover set

Item	Description
Safety information	Observe the "Safety Information" section in the first chapter, "General Information".
Crane suspension	A crane suspension for moving the gantry via crane is available and has to be ordered directly from the supplier. Order procurement: - Order designation: "Anschlagpunkt" - Order number: 166199 - Supplier/order address: PFEIFER Seil- und Hebetechnik GmbH Dr.-Karl-Lenz-Strasse 66 87700 Memmingen Germany NOTES: - The eye hook must be ordered with the German designation. - For safety reasons, we recommend using safety cords located at the transport device between the crane hook and the eyelet. It must not be used for craning! The project manager has to provide a torque wrench to tighten the hexagon nut (item 7) with 250 Nm. Fig. 74: Crane suspension for the gantry (schematic illustration) (A) Front view (B) Dimensions of the eye hook (1) Location of the eye hook (barycenter of the gantry) (2) Loops of the safety cord; attached to the crane hook (3) Recommended safety cords, length approx. 3 m; on-site (4) Eyelet on the transport device (5) Transport device (6) Thread size (metric) (7) Hexagonal nut, size 36; torque 250 Nm

Adaptor kit for the transport device:

Usually, a transport device is used with more than one CT system. For certain CT systems, special adaptor kits are required for attaching the gantries to the commonly used transport devices.

- For new CT systems, adaptors are always included in the scope of delivery and are preinstalled on the gantry where required.
- The system-specific adaptor kit needs to be ordered together with the transport device (for relocations, dismantling, etc.).

Order procurement:

The Transport device must be ordered separately for all non-European countries.

Functions:

- The transport height is adjustable in both directions. Depending on the type of transport device, the gantry can be lifted up and down via a spindle or using a hydraulic lifting system.
- The transport rollers can be swiveled in and out.

Function	Use
Swiveling out	During standard transport
Swiveling in	The transport width can be reduced when moving through narrow passages during transport.

CAUTION

[hz_pgdoc_F13G01U16M02]

Transport with the rollers swivelled in is permissible only when narrow places make it impossible to transport the system with the rollers swivelled out.

Tipping hazard!

» As soon as the system has passed through narrow passages, the transport rollers must be swivelled out again.

hz_pgdoc_F13G01U16M02

Transport dimensions:

The total gantry transport dimensions are defined by the overall dimensions of the gantry with the attached transport device.

Transport weights:


The total transport weight always includes the weight of the transport device with the adaptor kit (if one is used).

Designation	Weight
Gantry with the transporting device (without packaging)	<2900 kg

Transporting device	280 kg
---------------------	--------

In special cases, it may be necessary to reduce the total transport weight. This can be done by removing parts of the cover. In this case, a service engineer is required to disassemble and reassemble the device.

Cover designation	Weight	Remarks
Front ring upper		Do not remove the upper front ring! Reason: The patient positioning lasers are integrated in the upper front ring. After removal of this cover, a costly and time-consuming readjustment of the lasers by factory support is required. NOTE: Prior agreement from CR-factory is unavoidable. Contact CR SCM Installation Management.
Sealing ring	1,00 kg	n.a.
Ring front bottom	5,70 kg	
Cover right	15,764 kg	
Side cover right	10,952 kg	
Cover front right	8,210 kg	
Arc front left	4,570 kg	
Arc front center	0,715 kg	
Baseboard front right	2,760 kg	
Connector Cover	0,172 kg	
Cover front left	3,792 kg	
Baseboard front center	1,760 kg	
Ring rear with GPA	24,40 kg	
Arc rear left	3,850 kg	
Cover left	12,707 kg	
Funnel	13,364 kg	
Side cover left	8,672 kg	
Baseboard rear right	2,630 kg	
Service cover rear	0,495 kg	
Cover rear left	4,432 kg	
Cover load hook	0,360 kg	
Caps Set Cover	0,05 kg	
Flap mechanism frontcover	7,60 kg	
GPA left	2,150 kg	
GPA right	2,10 kg	

Pos.	Component
1	Ring front top
2	Sealing ring
3	Ring front bottom
4	Cover right
5	Side cover right
6	Cover front right
7	Arc front left
8	Arc front center
9	Baseboard front right
10	Connector Cover
11	Cover front left
12	Baseboard front center
13	Cover rear right
14	Ring rear with GPA
15	Arc rear left
16	Cover left
17	Funnel
18	Side cover left
19	Baseboard rear right
20	Service cover rear
21	Cover rear left
22	Cover load hook
24	Caps Set Cover
25	Flap mechanism frontcover
27	GPA left
28	GPA right

Floor load during transport:

Design access floors for a weight capacity of min. 400 kg per slab/plate. During gantry transport, the load may be higher at certain individual points (3-point load, e.g. due to uneven flooring).

Avoid damage to the on-site flooring. If necessary, cover the transport route with metal sheets to distribute the load and protect the customer's flooring.

15.5. TRANSPORT

Standard transport:

Standard transport means that the rollers of the transport device are swiveled out.

Transport through narrow passages:

The rollers of the transport device can be swiveled in for transporting items through narrow spaces, such as narrow doors.

CAUTION

[hz_pgdoc_F13G01U16M02]

Transport with the rollers swivelled in is permissible only when narrow places make it impossible to transport the system with the rollers swivelled out.

Tipping hazard!

- » As soon as the system has passed through narrow passages, the transport rollers must be swivelled out again.

hz_pgdoc_F13G01U16M02

Crane suspension for the gantry:

Using the transport device as a crane is not permitted!

If craning is necessary, a crane must be provided and set up on-site in accordance with the applicable local and national regulations in your country.

Contract a structural engineer if necessary.

15.6. MOVING THE PATIENT TABLE

The patient tables of CT systems are fully assembled and delivered on a wooden pallet. A ramp to move the PHS from the pallet is included in the delivery.

A transport device (transport device with rollers) is premounted at the table base. It must be removed after installation and stored at the customer site. The transport route must be smooth and level.

Store the transport device at the customer site for possible subsequent use by service or moves within the hospital, and so on.

15.7. MOVING THE IRS TOWERS AND IRS CABINETS

The IRSmx towers and IRS cabinets are delivered on a transport pallet to facilitate their transport through buildings.
With integrated rollers at the housing bottom, the IRSmx towers and IRS cabinets can be moved to their final installation location within the equipment room.

16. Storage and service life.

Climatic conditions for storage:

Designation	Range
Temperature	-20° C to +50° C
Relative humidity	10% to 90%
Temperature gradient	10 K/h
Barometric pressure	700 hPa to 1060 hPa

Service life:

Service life of the medical device is 10 years.

17. Procedure of recycling and disposal.

 **CAUTION**

Improper disposal of the system or parts of the system!
Pollution of the environment.

- ◆ System components hazardous to persons or the environment must be disposed of with care and in compliance with legally binding ordinances.
- ◆ Examples of environmentally relevant components are accumulators and batteries, transformers, capacitors, monitor picture tubes, and phantoms.
- ◆ For details contact your local customer service representative or your Siemens regional office.

17.1. GENERAL INFORMATION

Handling of hazardous materials:

Current recycling and waste management regulations assign product responsibility to the manufacturer/seller. Compliance with these product regulations includes evaluation of recyclables and environmentally safe disposal of materials and components.

This document provides information about the type, quantity and location of hazardous materials in the as well as instructions how to remove the relevant components. The disposal instructions have to be adhered to when disposing of the system. When disposing of the system and/or parts, please note that repairable parts have to be returned for reuse.

 **CAUTION**

[hz_serdoc_F14G08U01M02]

Hazardous materials have to be properly disposed of by an authorized waste management company.

Risk of environmental pollution.

- » Dispose of hazardous materials in accordance with local regulations.

hz_serdoc_F14G08U01M02

17.2. LIST OF HAZARDOUS MATERIALS

Components	Hazardous materials	Quantity
------------	---------------------	----------

Gantry:	Lead	15640 grams
	In tube collimators:	
	- X-ray tube assemblies	2 x 7400 grams
	- Tube collimators	2 x 420 grams
	Insulating oil	56,8 liters (total)
	Contained in:	
Patient table (standard) or Multi-purpose patient table	- High voltage transformers	2 x 17,1 liters
	- X-ray tube assemblies	2 x 4,5 liters
	- Tube cooling units	2 x 6,8 liters
	Grease in main bearing	150 grams
	Antimony trioxide	Traces ¹⁾
	Arsenic	
	Grease	
	In spindle of vertical drive, upper and lower guide rails of tabletop and top support	Small amount
Power distribution cabinet includes: Uninterruptible power supply (UPS): (located inside PDC)	Rechargeable batteries In UPS	1 battery pack
Water-air split cooling system (WCS-w-a-split): (Indoor unit) (Outdoor unit)	Coolant liquid (R134a) (in refrigerant circuit of the indoor unit)	2600 grams
	Compressor oil (in compressor of indoor unit)	2900 grams
	Antifreeze (in cooling water circuit between indoor unit and outdoor unit)	40% Glycol Amount depends on actual water hose length
Image control system (ICS-tower):	Lithium battery (on motherboard)	1 unit
Image evaluation system (IES-tower):	Lithium battery (on motherboard)	1 unit
Image reconstruction system (IRS):	Lithium battery (on motherboard)	1 unit
Color monitor:	Antimony trioxide / Arsen in LCD display	Traces ¹⁾

Interventional Module (OPTION): includes: Interventional Panel (IVP):	Rechargeable batteries in IVP	4 units (AA-type)
Keyboard & Mouse	Not applicable	Not applicable
Control box:	Not applicable	Not applicable
¹⁾ small amounts in the gantry covers, control panels and LCD-display		

17.3. DISPOSAL OF HAZARDOUS MATERIALS

General information:

WARNING

[hz_serdoc_F13G01U12M01]

Failure to observe the guidelines specified in this document may lead to serious injury or loss of life.

Risk of accident and/or injury!

- » The guidelines contained in this document have to be strictly followed.

hz_serdoc_F13G01U12M01

WARNING

[hz_serdoc_F13G01U05M03]

Dangerous voltages are present when the system is switched on. Dangerous voltages may be present even when the system is switched off due to capacity power.

Risk of electric shock!

- » Observe power-off instructions and discharge wait times (allow at least five minutes discharge time after the last scan for all involved HV and UDC/UPS parts).
- » Secure the system against unintended switch-on (e.g. block breaker and/or mark breaker against switch on).
- » Ensure via measurements that all voltages are switched off.
- » Connect the stationary and rotating parts of the gantry to a protective conductor prior to working in the gantry.

hz_serdoc_F13G01U05M03

WARNING

[hz_serdoc_F12G02U07M02]

The UPS contains its own energy source (batteries). The parts supplied with UPS power may still have power even when the UPS is not connected to an AC supply (power off status).

Risk of electrical shock from the UPS!

- » Switch off UPS power by switching off the relevant circuit breakers.
- » Secure the system against unintended switch-on (e.g. block braker against switching on and/or mark breaker against switch on).

hz_serdoc_F12G02U07M02

CAUTION

[hz_serdoc_F13G01U01M08]

Unexpected gantry rotation during replacement of parts.

Risk of accident and injury!

- » Secure the gantry with the "Rotation safety bolt" before replacing any parts in the rotating part of the gantry.
- » Do not remove the "Rotation safety bolt" until all parts are reinstalled and secured.

hz_serdoc_F13G01U01M08

CAUTION

[hz_serdoc_F13G07U01M01]

Handling parts of the system that may have come into contact with patients may lead to infection caused by blood-borne pathogens.

Infection caused by blood-borne pathogens!

- » Take appropriate precautions against exposure to blood-borne pathogens (e.g. wear gloves).

hz_serdoc_F13G07U01M01

CAUTION

[hz_serdoc_F14G08U01M02]

Hazardous materials have to be properly disposed of by an authorized waste management company.

Risk of environmental pollution.

» Dispose of hazardous materials in accordance with local regulations.

hz_serdoc_F14G08U01M02

18. Maintenance.

On request, you can have maintenance work performed regularly by the Siemens Customer Service. Please contact your Siemens Service to agree on a maintenance contract for your system.

Routine checks:

WHAT:	WHEN:
Image quality test	Daily
Camera test	Daily
Function checks	Daily
Constancy test*	Monthly

*We recommend that the constancy test be performed by Siemens service personnel in conjunction with a maintenance contract.

Periodic maintenance activities are subdivided into:

WHAT:	WHEN:
Partial maintenance	Every 6 months
Full maintenance	Every 12 months
Preventive replacement of safety-relevant parts	Every 2 years

This work must be performed by qualified and authorized service engineers. In this context, qualified means that the service engineer has been trained and has had practical experience in the required routines so that he/she is capable of performing service work on the system. Authorized means that the service engineer is recognized by the system owner as a qualified service engineer and is therefore authorized to perform maintenance work on the system.

Upon system start-up, we recommend designating one employee who is responsible for ensuring that routine checks and preventive inspection and maintenance work are performed. The same employee should also be responsible for archiving all protocols in the "System Owner Manual. In addition to repair services, Siemens also offers its customers complete maintenance service, including inspection and preventive maintenance on Siemens systems. These services can be provided either on a service call basis or in the form of a flexible maintenance contract. If you have not received a quotation for the performance of preventive maintenance and inspection from our Siemens UPTIME Services, your local Siemens sales representative will gladly provide this information for you.

For service and maintenance on the territory of Russian Federation please use the following contacts:

*Siemens Healthcare LLC
Dubinskaja str. 96
115093 Moscow, Russia
Tel.: +7 495 737 10 00*

19. Calibration.

System related artifacts:

System related artifacts may be due to the following:

- Uncalibrated system;
- Detector deviations.

 **CAUTION**

When you restart the system, the detector has not yet reached operating temperature!

Wrong diagnosis due to image artifacts.

- ◆ Calibrate the system as part of the checkup. Repeat calibration (**Setup > Calibration**) if ring artifacts occur.

Calibration:

If the CT system is not properly set up (for example, no checkup performed) CT scale displacements and inhomogeneities can result. This may result just after switching on the unit, during warm-up to operational temperature (calibration), or due to an extended length of service of the X-ray tube. The same applies to defects in the measurement system. The daily quality checkup usually detects and reports such discrepancies. Other system defects or calibration deviations cause streaks or (partial) rings to appear in the CT image.

Detector deviations:

In CT systems with a rotating combined tube-detector system even minimal deviations of individual detector channels from the original calibration level may lead to rings or partial ring artifact structures in the CT image. The closer such channels lie to the detector center, the greater are these effects.

CAUTION

When you restart the system, the detector has not yet reached operating temperature!

Wrong diagnosis due to image artifacts.

- ◆ Calibrate the system as part of the checkup. Repeat calibration (**Setup > Calibration**) if ring artifacts occur.

In an extreme case rings may appear concentrated as a blurred spot in the center of the scan field (the scan field center appears at the image center if the **Center X** and **Y** image parameters are both zero). Such rings are easily recognizable, but a blurred spot in the center of the scan field might lead to inaccurate diagnoses.

20. Warranty.

1. We warrant that the deliveries are free from material defects and lack none of the warranted characteristics at the time of risk transfer. Warranted characteristics are only those characteristics explicitly designated as such in the specifications. Our relevant Special Supply and Warranty Conditions shall apply in respect of products for which varying warranty regulations are agreed for product-specific reasons, such as high-vacuum elements or flat detectors.
2. The warranty period is 12 months from installation and commissioning, but shall expire no later than 15 months after delivery in accordance with the agreed Incoterms. The period of warranty for defects in respect of spare parts purchased from us is 6 months from the date of risk transfer.
3. Warranty claims shall not arise from trivial deviation from the agreed state or condition, in the case of non-reproducible software errors or in the course of normal wear and tear. Neither shall they arise from damages resulting from improper or negligent handling, assembly or commissioning by the Purchaser or third parties, excessive demand, improper use, failure to comply with operating instructions, use of inappropriate operating materials, incorrect structural prerequisites or technical specifications of the Purchaser or third parties, or due to special external influences not provided for in the contract. If modifications or repairs are performed by the Purchaser or a third party, no warranty claims shall exist in respect of such modifications or repairs or the consequences arising therefrom. Unless otherwise agreed, we accept no liability for material defects in used products.
4. Immediately after receiving the deliveries, the Purchaser shall examine them from defects. Obvious defects shall be reported immediately after shipment or installation or at the latest within fourteen (14) calendar days of receipt of the deliveries. Hidden defects shall be reported as soon as they are identified. Complaints shall be raised in writing, specifying the order data, the invoice and shipment numbers, and the identification number of the defective deliveries. If the Purchaser fails to provide notification in the correct form and in good time, the deliveries shall be deemed approved, subject to the reservation that hidden defects may be reported at a later time. The Purchaser shall not refuse to accept delivery due to trivial defects.

5. Purchaser shall provide Siemens with both on-site and remote access to the Products to the extent legally permitted and technically feasible. The remote access shall be provided through Purchaser's Internet Connection as is reasonably necessary for Siemens to provide services under this Contract, such as rectification of defects. The Purchaser agrees to allow connection to Siemens' remote service system. In the case of remote services, Siemens shall be entitled to access the data available with the Purchaser for the Products serviced by remote service.

6. Deliveries that exhibit material defects or lack an explicitly warranted characteristic at the time of risk transfer shall, at our discretion and within the normal processing times, be either repaired free of charge or replaced with non-defective deliveries. We shall assume ownership of any defective parts replaced. In the event of defects in the software supplied to Purchaser in connection with the Products – substantial non-conformance with the relevant software specification- Siemens provides Purchaser with an updated version of the software product eliminating such defect, which is or becomes available to Siemens. In the event modified standard software or developed individual software is defective, Siemens will until the provision of an updated version of the software eliminating the defect provide Purchaser with an acceptable by-pass or error eliminating solution, if this is possible with reasonable expenditure and if otherwise the commercial use of the Product is substantially affected.

7. Unless otherwise expressly agreed, we shall deliver repaired or replacement parts at our discretion FOB German port or FCA German airport in accordance with incoterms 2010. We shall determine the port or the German airport to be used. We accept no responsibility for sending back or storing deliveries returned to us without our prior consent. For quality reasons, goods may be returned or taken back only in agreement with our order processing department.

8. If repair or replacement delivery is unsuccessful or if we allow a reasonable extension period to expire without rectifying the defect, the Purchaser may demand a price reduction or, if the defective delivery or the part of the delivery that is free from defects is of no use to the Purchaser, may withdraw from the contract.

9. No new warranty period shall commence as a result of repair, replacement or substitute performance/delivery.

10. When defects are reported, payments by the Purchaser may be withheld to an extent reasonably proportionate to the defects identified. However, the Purchaser may withhold payment only if a justifiable complaint has been brought. If the complaint made is unfounded, we are entitled to demand reimbursement from the Purchaser of our costs incurred thereby.

11. Further claims and rights of the Purchaser beyond those provided in this Item 10 asserted against us due to material defects or the lack of warranted features are excluded.

21. Contacts for reclamation.

In case of any questions or problems with handling the system, please use the following contacts:

Contacts for questions and reclamations in Russia:

Siemens Healthcare LLC

Dubiniskaya street 96

115093 Moscow,

Russia

Tel.: +7 495 737 12 52

Fax: +7 495 737 13 20

Customer support center

8 800 550 50 22

(24 hour free line)

service.ru@siemens-healthineers.com

Contacts for questions and reclamations in the country of Legal manufacturer:

Siemens Healthcare GmbH

Henkestr. 127,

91052 Erlangen

Germany

Tel.: +49 9131 84-0

This document has been printed in Germany.

Siemens Healthcare GmbH

Rüdiger Körner
Head of Regulatory Affairs
Business Line Computed Tomography


Deed No. N 1271 /2021 Th

I, Marcel Neumair, notary in Forchheim, notarize that the document overleaf was signed in my presence by

Mr. Ralf Rüdiger Körner,
born on September 24th, 1957,
with business address in D-91301 Forchheim, Siemensstraße 1,
personally known to me, notary,

acting according to his declaration for
Siemens Healthcare GmbH.

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