Medtronic

18000 Devonshire Street
Northridge, CA 91325 USA
www.medtronic.com
Tel 818-576-5211

ATTESTATION STATEMENT

As a responsible member of Medtronic MiniMed, I hereby attest that I believe the Instructions for Use attached herein, to be a true and accurate copy.

Katherine Rojas
Regulatory Affairs Specialist

Katherine Rojas July 16, 2021
Signature

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document


State of California
County of Ventura

On July 16th 2021 before me, Kusum Manwani
(date) (Insert name and title of the officer)

personally appeared Katherine Rojas, who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the same in his/her/their authorized capacity(ies); and that by his/her/their signature(s) on the instrument the person(s), or the entity upon behalf of which the person(s) acted, executed the instrument. I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

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Medtronic

Contacts:

Africa: Medtronic South Africa and Southern Africa

Office Reception Tel: +27(0) 11 260 9300

Diabetes: 24/7 Helpline: 0800 633 7867

Sub-Sahara 24/7 Helpline: +27(0) 11 260 9490

Albania: Net Electronics Albania

Tel: +355 697070121

Argentina: Corpomedica S.A.

Tel: +(11) 4 814 1333

Medtronic Directo 24/7:

+0800 333 0752

Armenia: Exiol LLC

Tel: +374 98 92 00 11

or +374 94 38 38 52

Australia: Medtronic Australasia Pty. Ltd.

Tel: 1800 668 670

Bangladesh: Sonargaon Healthcare Pvt Ltd.

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Malaysia: Medtronic International Ltd.
Tel: +603 7946 9000

Middle East and North Africa: Regional Office
Tel: +961-1-370 670

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Sri Lanka: Swiss Biogenics Ltd.
Mobile: (+91)-9003077499
or (+94)-777256760

Suomi: Medtronic Finland Oy
Tel: +358 20 7281 200
Help line: +358 400 100 313

Sverige: Medtronic AB
Tel: +46 8 568 585 20
Fax: +46 8 568 585 11

Taiwan: Medtronic (Taiwan) Ltd.
Tel: 02-21836000
Toll free: +886-800-005285

Thailand: Medtronic (Thailand) Ltd.
Tel: +662 232 7400

Türkiye: Medtronic Medikal Teknoloji
Ticaret Ltd. Sirketi.
Tel: +90 216 4694330

USA: Medtronic Diabetes Global Headquarters
24-Hour Technical Support: +1-800-646-4633
To order supplies: +1-800-843-6687

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
The Guardian™ Link (3) transmitter with Bluetooth™* wireless technology is a component of the continuous glucose monitoring (CGM) system for the MiniMed™ insulin pump systems with smart device connectivity.

Guardian™ Link (3) transmitter kit (MMT-7910) components

A complete transmitter kit includes the following components:

- Guardian™ Link (3) transmitter (MMT-7911)
- Two testers (MMT-7736L)
- Charger (MMT-7715)
- One-press serter (MMT-7512)

Intended purpose of the device

The Guardian™ Link (3) transmitter (MMT-7911) is a rechargeable device and powers the glucose sensor, collects and calculates the sensor data, and sends the data to a compatible MiniMed™ insulin pump system with smart device connectivity for the management of diabetes mellitus. The transmitter is only compatible with the Guardian™ Sensor (3) glucose sensor (MMT-7020) and is indicated for single-patient or multiple-patient use.

Contraindications

None known.

Warnings

- Do not use the transmitter adjacent to other electrical equipment that may cause interference with the normal system operation. Other electrical equipment that may compromise normal system operation has been contraindicated. For more information on electrical equipment that may compromise normal system operation, see *Exposure to magnetic fields and radiation*, on page 2.


- Always refer to the sensor user guide for all precautions, warnings, and instructions relating to the sensor. Not referring to the sensor user guide can result in serious injury or damage to the sensor.
- Do not allow children to put small parts in their mouth. This product poses a choking hazard for young children.
- Do not change or modify the device unless expressly approved by Medtronic Diabetes. Modifying the device can cause serious injury, interfere with your ability to operate the device, and void your warranty.
- Do not use the tester if it comes in contact with blood. Touching blood can cause infection. Dispose of the tester according to the local regulations for medical waste disposal, or contact your healthcare professional for disposal information.
- Bleeding may occur after inserting the sensor. Always make sure that the site is not bleeding before connecting the transmitter to the sensor. Blood can get into the transmitter connector and damage the device. Discard the device if damaged. If bleeding occurs, apply steady pressure with a sterile gauze or clean cloth at the insertion site until bleeding stops. After bleeding stops, connect the transmitter to the sensor.
- Contact your local representative if you experience any adverse reactions associated with the transmitter or sensor. Adverse reactions can cause serious injury.
- Do not discard the transmitter in a medical waste container or expose it to extreme heat. The transmitter contains a battery that may ignite and result in serious injury.

Exposure to magnetic fields and radiation

- Do not expose your transmitter to Magnetic Resonance Imaging (MRI) equipment, diathermy devices, or other devices that generate strong magnetic fields (for example, x-ray, CT scan, or other types of radiation). Exposure to a strong magnetic field has not been evaluated and can cause the device to malfunction, result in serious injury, or be unsafe. If your transmitter is exposed to a strong magnetic field, discontinue use and contact your local representative for further assistance.
- Always remove your sensor and transmitter before entering a room that has x-ray, MRI, diathermy, or CT scan equipment. Exposure to a strong magnetic field has not been evaluated and can cause the device to malfunction, result in serious injury, or be unsafe. If your sensor or transmitter is exposed to a strong magnetic field, discontinue use and contact your local representative for further assistance.
- Always carry the Medical emergency card provided with your device when you are traveling. The Medical emergency card provides critical information about airport security systems and using your transmitter safely on an airplane. Not following the guidance on the Medical emergency card could result in serious injury.


Precautions

- Do not attempt to use the Guardian™ Link (3) transmitter (MMT-7911) with a MiniMed™ insulin pump without smart device connectivity. Only a MiniMed™ insulin pump with smart device connectivity can communicate with the Guardian™ Link (3) transmitter (MMT-7911).
- Only use the Guardian™ Sensor (3) glucose sensor (MMT-7020) with the transmitter. Do not use any other sensor. Other sensors are not intended for use with the transmitter and will damage the transmitter and the sensor.
- Only use the green colored tester (MMT-7736L) with the transmitter. Pockets on the transmitter are visible when connected to the tester. Do not use any other test plug. Other test plugs are not intended for use with the transmitter and will damage the transmitter and the tester.

- Always use the tester when cleaning the transmitter. Do not use any other test plug with the transmitter. Use of another test plug can allow water to get into the transmitter or can prevent proper cleaning. Water can damage the transmitter.
- Do not twist the tester or sensor while attached to the transmitter. Twisting the tester or sensor will damage the transmitter.
- Do not allow the tester to come in contact with any liquid when not connected to the transmitter. A wet tester can damage the transmitter.
- Do not allow the transmitter to come in contact with any liquid when not connected to a sensor or to the tester. Moisture will damage the transmitter and a wet transmitter can damage the sensor.

- Do not clean the O-rings on the tester with any substances. Cleaning the O-rings can damage the tester.

IEC 60601-1-2:2014, 4th Edition; Special EMC Precautions for Medical Electrical Equipment

- 1 Special Precautions regarding Electromagnetic Compatibility (EMC): This body worn device is intended to be operated within a reasonable residential, domestic, public or work environment where common levels of radiated "E" (V/m) or "H" fields (A/m) exist, such as cellular phones, wireless technology, electric can openers, microwave and induction ovens. This device generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the provided instructions, may cause harmful interference to radio communications.
- 2 Portable and mobile RF communications equipment can affect medical electrical equipment. If you encounter RF interference from a mobile or stationary RF transmitter, move away from the RF transmitter that is causing the interference.
- 3 Be careful when you use your transmitter closer than 30 cm (12 in) to portable radio frequency (RF) equipment or electrical equipment. If you must use your transmitter next to portable RF equipment or electrical equipment, observe the transmitter to verify correct system operation. Degradation of the performance of the transmitter could result.

Assistance

Contact your local representative if you need a copy of a MiniMed™ system user guide.

Preparing your transmitter


The transmitter contains a non-replaceable, rechargeable battery that you can recharge as needed with the charger. The transmitter needs to be charged before you use it. The charger has a green light that shows the charging status and a red light that communicates any problems during charging. If you see a red light, see *Troubleshooting*, on page 19. The charger needs one AAA alkaline battery.

Note: If the battery is installed incorrectly or is low, the charger will not work. Repeat the battery installation steps using a new battery.

Installing a battery in the charger

To install a battery in the charger:

- 1 Push the battery cover in and slide it off (as shown in the image in step 3).
- 2 Insert a new AAA alkaline battery. Make sure the + and - symbols on the battery align with these same symbols shown on the charger.
- 3 Slide the cover back on the charger until it clicks into place.


Charging the transmitter

CAUTION: Always charge the transmitter before inserting your sensor. A depleted transmitter does not function. A fully charged transmitter works at least seven days without recharging. A depleted transmitter can take up to two hours to recharge.

CAUTION: Do not store the transmitter on the charger for more than 60 days. Disconnect and reconnect to the charger to re-charge again before use. If the transmitter is left on the charger for more than 60 days, the transmitter battery will be permanently damaged.

To charge the transmitter:

- 1 Push the transmitter and the charger together to connect the transmitter to the charger.
- 2 Within 10 seconds after the transmitter is connected, a green light on the charger will flash for one to two seconds as the charger powers on. For the rest of the charging time, the green light on the charger will continue to flash in a

pattern of four flashes with a pause between the four flashes.

- 3 When charging is complete, the green light on the charger will stay on, without flashing, for 15 to 20 seconds and then turn off.
- 4 After the green charger light turns off, disconnect the transmitter from the charger. The green light on the transmitter starts to flash.

Pairing your transmitter

Always refer to the system user guide for instructions on how to pair your transmitter with your pump. The pump and the transmitter must be paired before data from the sensor can be sent to the pump. The pump and the transmitter are only required to be paired once. There is no need to pair the pump with the transmitter again when you insert a new sensor.

Inserting the sensor


Always refer to your sensor user guide for instructions on how to insert the sensor.

Connecting the transmitter to the sensor

Before proceeding, have your MiniMed™ insulin pump system user guide available.

To connect the transmitter to the sensor:

- 1 After the sensor is inserted, consult your sensor user guide for details on how to apply the required tape before connecting the transmitter.
- 2 Hold the rounded end of the inserted sensor to prevent it from moving during connection.
- 3 Hold the transmitter as shown. Line up the two notches on the transmitter with the side arms of the sensor. The flat side of the transmitter should face the skin.
- 4 Slide the transmitter onto the sensor connector until the sensor arms snap into the notches on the transmitter. If the transmitter is properly connected, and if the sensor has had enough time to become hydrated with the interstitial fluid, the green light on the transmitter will flash 6 times.

Note: If the transmitter does not flash, see *Troubleshooting*, on page 19.


- 5 When the transmitter light flashes green after connecting to the sensor, use your pump to start the sensor. For more instructions, see your system user guide.
- 6 Attach the adhesive tab of the sensor to the transmitter.
- 7 After the transmitter is connected, consult your sensor user guide for details on applying the required tape.
- 8 Follow the instructions that appear on the pump screen or in your system user guide.

Disconnecting the transmitter from the sensor

Before proceeding, have your MiniMed™ insulin pump system user guide available.

To disconnect the transmitter from the sensor:

- 1 Carefully remove any tape from the transmitter and sensor.
- 2 Remove the adhesive tab from the top of the transmitter.
- 3 Hold the transmitter as shown, and pinch the flexible side arms of the sensor between your thumb and forefinger.
- 4 Gently pull the transmitter away from the sensor.
- 5 Follow the instructions that appear on the pump or in your system user guide.

Removing the sensor

Always refer to the sensor user guide for instructions on how to remove the sensor.

Reconnecting the transmitter to a sensor that is already inserted

You can reconnect your transmitter to the sensor you are currently using. Simply connect your transmitter to the sensor that is already inserted. When the pump detects the transmitter, confirm that you want to Reconnect Sensor. It may take a few seconds to establish a connection when reconnecting a sensor. Reattach the adhesive tab of the sensor to the transmitter and reapply any required tape. When you reconnect a sensor, the sensor will go through another warm-up period before you can calibrate it.


Tester

The tester is used to test the transmitter to make sure it is working. The tester is also used as a required component to create a waterproof seal when cleaning the transmitter. Properly connecting the tester to the transmitter ensures that fluids do not come in contact with the connector pins inside the transmitter. Fluids can cause connector pins to corrode and affect the performance of the transmitter.

Do not twist the tester while it is attached to the transmitter. This will damage the transmitter.


The tester can be used for one year. If you continue to use the tester for more than one year, the connector pins inside the transmitter could be damaged, because the tester cannot continue to provide a waterproof seal. For instructions on how to check the connector pins, see *Inspecting the transmitter connector pins*, on page 8.

CAUTION: Only use the green colored tester (MMT-7736L) with the transmitter. Pockets on the transmitter are visible when connected to the tester. Do not use any other test plug. Other test plugs are not intended for use with the transmitter and will damage the transmitter and the tester.

Inspecting the transmitter connector pins

This image is an example of how the connector pins should look.

Look inside the connector opening of the transmitter to make sure that the connector pins are not damaged or corroded. If the connector pins are damaged or corroded, the transmitter cannot communicate with the charger or the pump. Contact your local representative. It may be time to replace your transmitter.

Also look for moisture inside the connector opening. If you see any moisture, allow the transmitter to dry for at least one hour. Moisture inside the connector opening could cause the transmitter to not work properly and could cause corrosion and damage over time.

Connecting the tester for testing or cleaning


Before proceeding, have your MiniMed™ insulin pump system user guide available.

To connect the tester:

- 1 Hold the transmitter and the tester as shown. Line up the flat side of the tester with the flat side of the transmitter.
- 2 Push the tester into the transmitter until the flexible side arms of the tester click into the notches on both sides of the transmitter.

When properly connected, the green light on the transmitter flashes 6 times.

- 3 To test the transmitter, check the sensor icon on the pump to ensure that the transmitter is sending a signal (see your system user guide).
- 4 To clean the transmitter, see *Cleaning the transmitter*, on page 9.
- 5 After testing or cleaning, disconnect the tester from the transmitter.


Disconnecting the tester

To disconnect the tester:

- 1 Hold the transmitter body as shown and pinch the side arms of the tester.
- 2 With the tester arms pinched, gently pull the transmitter away from the tester.

Note: To save transmitter battery life, do NOT leave the tester connected after cleaning or testing.

Cleaning the transmitter

The transmitter is intended for personal use at home (single-patient use) or for use in healthcare facilities (multiple-patient use). Single-patient use requires cleaning after each use, while multiple-patient use requires cleaning and disinfection after each use. When using the transmitter in a healthcare facility, always follow the cleaning and disinfecting procedure for multiple-patient use.

WARNING: Do not discard the transmitter in a medical waste container or expose it to extreme heat. The transmitter contains a battery that may ignite and result in serious injury.

Note: The tester is a required component for cleaning the transmitter. For details, see *Tester*, on page 7.

CAUTION: Do not use automated washer-disinfector to clean or disinfect the device. Using automated washer-disinfector to clean or disinfect the device will cause damage to the transmitter.

For single-patient use

Always clean the transmitter after each use.

To clean the transmitter, you need the following materials:

- mild liquid soap
- soft-bristled toddler toothbrush
- container
- clean, lint-free dry cloths


Use life

The transmitter can be cleaned up to 122 times or for one year, whichever comes first. Discard the transmitter at this point. If you continue to use the transmitter beyond 122 times or one year, the cleaning process may damage the device. Contact your local representative to order a new transmitter.


WARNING: Do not use the device if you see any cracking, flaking, or damage to the housing. Cracking, flaking, or damage to the housing are signs of deterioration. Deterioration of the housing can affect the ability to properly clean the transmitter and result in serious injury. Call your local representative and discard the device according to local regulations for battery disposal (non-incineration), or contact your healthcare professional for disposal information.

To clean the transmitter:


- 1 Wash your hands thoroughly.
- 2 Attach the tester to the transmitter to create a waterproof seal.

- 3 If there is adhesive residue on the transmitter, see *Removing adhesive residue, on page 17*.
- 4 Rinse the transmitter under room temperature tap water for at least one minute, and until visibly clean. Make sure all hard-to-reach areas are rinsed completely.


- 5 Prepare a mild liquid soap solution using 5 milliliters (1 teaspoon) of mild liquid soap per 3.8 liters (1 gallon) of room temperature tap water.
- 6 With the tester still attached, submerge the transmitter in the mild liquid soap solution and soak for one minute.


- 7 Holding the tester, brush the entire surface of the transmitter using a soft-bristled toddler toothbrush. Make sure to brush all hard-to-reach areas until visibly clean.

- 8 Rinse the transmitter under running room temperature tap water for at least one minute, and until all visible liquid soap is gone.**


- 9 Dry the transmitter and tester with a clean, dry cloth.**

- 10 Place the transmitter and tester on a clean, dry cloth and air dry them completely.**

- 11 Disconnect the tester from the transmitter by gently squeezing the arms of the tester.**

For multiple-patient use

When using the transmitter in a healthcare facility, always clean and disinfect the transmitter after each use.

WARNING: You must adhere to Standard Precautions when handling or using this device. All parts of the system should be considered potentially infectious and are capable of transmitting blood-borne pathogens between patients and healthcare professionals.

The transmitter must be disinfected after use on each patient. This system may only be used for testing multiple patients when Standard Precautions and disinfection procedures provided by Medtronic Diabetes are followed.

To clean the transmitter, you need the following materials:

- gloves
- mild liquid soap
- soft-bristled toddler toothbrush
- 8.25% bleach
- two containers
- clean and lint-free dry cloths

Use life


The transmitter can be cleaned and disinfected up to 122 times or for one year, whichever comes first. Discard the transmitter at this point. If you continue to use the transmitter beyond 122 times or one year, the cleaning and disinfection process may damage the device. Contact Medtronic to order a new transmitter.

To clean and disinfect the transmitter:


- 1 Wash your hands and put on gloves.
- 2 Inspect the inside of the transmitter's connector opening for any sign of body fluid. For instructions on how to inspect the connector pins, see *Inspecting the transmitter connector pins*, on page 8.

CAUTION: The person inspecting the transmitter must have sufficient vision that enables him or her to see small drops of body fluid or debris.


WARNING: If you see any body fluid in the connector opening, you must discard the transmitter. Because the transmitter contains a battery, do not discard in a bio-waste container. Instead, continue to clean and disinfect the transmitter, and then discard according to local regulations for battery disposal (non-incineration).


3 Attach the tester to the transmitter to create a waterproof seal.

- 4 If there is adhesive residue on the transmitter, see *Removing adhesive residue*, on page 17.**
- 5 Rinse the transmitter under room temperature tap water for at least one minute and until visibly clean. Make sure all hard-to-reach areas are rinsed completely.**


- 6 Prepare a mild liquid soap solution using 5 milliliters (1 teaspoon) of mild liquid soap per 3.8 liters (1 gallon) of room temperature tap water. Make sure to prepare a fresh solution for each use.
- 7 With the tester still attached, submerge the transmitter in the mild liquid soap solution and soak for one minute.


- 8 Holding the tester, brush the entire surface of the transmitter using a soft-bristled toddler toothbrush. Make sure to brush all hard-to-reach areas until visibly clean.

- 9 Rinse the transmitter under running room temperature tap water for at least one minute, and until all visible liquid soap is gone.


10 Dry the transmitter and tester with a clean, dry cloth.

- 11 Prepare a 1:10 bleach solution by using one (1) part 8.25% bleach to nine (9) parts water, for a final concentration of 0.8%. Make sure to prepare a fresh solution for each use.**
- 12 Make sure that you have completed the previous cleaning steps before disinfection. With the tester still attached, soak the transmitter in the bleach solution for 20 minutes.**


13 Rinse the transmitter under running room temperature tap water for three minutes.

14 Place the transmitter and tester on a clean, dry cloth and air dry them completely.

WARNING: If you saw any body fluid inside the connector opening on earlier inspection, you must now discard the transmitter with tester still attached, according to local regulations for battery disposal (non-incineration).

- 15 Disconnect the tester from the transmitter by gently squeezing the arms of the tester.

- 16 Inspect the housing of the transmitter for any signs of cracking, flaking, or damage. If you see any of these signs, you must now discard the disinfected transmitter according to local regulations for battery disposal (non-incineration).

WARNING: Do not use the device if you see any cracking, flaking, or damage to the housing. Cracking, flaking, or damage to the housing are signs of deterioration. Deterioration of the housing can affect the ability to properly clean the transmitter and result in serious injury. Call your local representative and discard the device according to local regulations for battery disposal (non-incineration), or contact your healthcare professional for disposal information.

- 17 Discard the used gloves and thoroughly wash hands with soap and water.

Removing adhesive residue


You may need to perform this procedure if there is adhesive residue present on the transmitter. If you visually inspect the transmitter and see adhesive residue on it, follow these instructions.

To remove adhesive residue, you need cotton swabs and a medical adhesive remover such as Detachol™*, which is a mineral spirit.

Note: During testing, Medtronic MiniMed used Detachol™* to remove the adhesive residue from the transmitter. Detachol™* is recommended for use but may not be available in all countries.

To remove adhesive residue:

- 1 Make sure the tester is attached to the transmitter.
- 2 Soak a cotton swab in the medical adhesive remover.
- 3 Hold the tester and gently rub the adhesive remover on the transmitter until the residue is removed.

- 4 Continue with the cleaning procedure. See *Cleaning the transmitter*, on page 9 for details.

Bathing and swimming

After the transmitter and sensor are connected, they form a waterproof seal to a depth of 2.4 meters (8 feet) for up to 30 minutes. You can shower and swim without removing them.

Cleaning the charger

This procedure is for general cleaning as required, based on physical appearance.

CAUTION: Do not immerse the charger in water or any other cleaning agent. The charger is not waterproof. Water can damage the charger and cause the device to malfunction.

WARNING: Dispose the charger according to the local regulations for battery disposal, or contact your healthcare professional for disposal information. The charger may ignite upon incineration.

To clean the charger:

- 1 Wash your hands thoroughly.

- 2 Use a damp cloth with mild cleaning solution, such as a dishwashing detergent, to clean any dirt or foreign material from the outside of the charger. Never use organic solvents, such as paint thinner or acetone, to clean the charger.
- 3 Place the charger on a clean, dry cloth and air dry for two to three minutes.

Troubleshooting

The following table contains troubleshooting information for the transmitter, charger, and tester. For more information about troubleshooting, see your system user guide.

Problem	Likely Cause(s)	Resolution
You connected the transmitter to the charger and no lights came on.	The transmitter connector pins are damaged or corroded. Your charger battery has no power or no battery is inserted.	1 Check the transmitter connector pins for damage or corrosion. For more information about your connector pins, see <i>Inspecting the transmitter connector pins</i>, on page 8. If the pins are damaged or corroded, contact your local representative. It may be time to replace your transmitter. 2 If there is no damage to the connector pins, replace the battery in the charger. For instructions on replacing your charger battery, see <i>Installing a battery in the charger</i>, on page 5.
During charging, the flashing green light on the charger turns off and you see a longer flashing red light on the charger.	Your charger battery is low on power.	Replace the battery in the charger. For instructions on replacing your charger battery, see <i>Installing a battery in the charger</i> , on page 5.
During charging, the flashing green light on the charger turns off and you see a series of quick flashing red lights on the charger for two seconds at a time.	Your transmitter is low on power.	1 Charge the transmitter continuously for one hour. If flashing does not stop, proceed to step 2. 2 Charge the transmitter continuously for eight hours. If flashing does not stop, call your local representative. It may be time to replace your transmitter.

Problem	Likely Cause(s)	Resolution
During charging, a mix of quick and long flashing red lights appear on the charger.	Your charger and your transmitter are low on power.	1 Replace the battery in the charger. For instructions on replacing your charger battery, see <i>Installing a battery in the charger</i>, on page 5. 2 Charge the transmitter continuously for one hour. If the quick flashing red lights do not stop, proceed to step 3. 3 Charge the transmitter continuously for eight hours. If flashing does not stop, call your local representative. It may be time to replace your transmitter.
The green light on the transmitter does not flash when you connect it to the sensor.	Your transmitter is not fully connected. Your transmitter is low on power. Your sensor is not properly inserted into your body.	1 Disconnect the transmitter from the sensor. 2 Wait for five seconds and reconnect them. If the green light still does not flash, proceed to step 3. 3 Fully charge the transmitter and connect it to the tester. If the green light still does not flash, see troubleshooting on "The green light on the transmitter does not flash when you connect it to the tester". If the green light flashes, proceed to step 4. 4 Disconnect the transmitter from the tester, wait at least five seconds, and connect the transmitter to the sensor. If the green light still does not flash, proceed to step 5. 5 The sensor may not be properly inserted into your body. Remove the sensor from your body and insert a new sensor.
The green light on the transmitter does not flash when you connect it to the tester.	Your transmitter is not fully connected. Your transmitter is low on power.	1 Check the connection between the transmitter and the tester. If the green light still does not flash, proceed to step 2. 2 Fully charge the transmitter. 3 Test the transmitter with the tester again. If you still do not see the green light flash, call your local representative. It may be time to replace your transmitter.

Problem	Likely Cause(s)	Resolution
Your transmitter battery does not last for seven days.	Your transmitter is not fully charged when you connect it to the sensor. The transmitter and pump frequently lose wireless connection.	1 Fully charge the transmitter before connecting it to the sensor. If the transmitter battery still does not last for the duration of one sensor use, proceed to step 2. 2 Move away from any equipment that can cause RF interference. For more information on RF interference, see the Radio Compliance Information sheet included with your pump. 3 Make sure your pump and your transmitter are located on the same side of your body to minimize any RF interference. If your fully charged transmitter battery continues to lose power before a full seven days, call your local representative. It may be time to replace your transmitter.
Your transmitter has lost connection with your pump.	Your pump is out of range. There is RF interference from other devices.	1 Move away from any equipment that can cause RF interference. For more information on RF interference, see the Radio Compliance information sheet included with your transmitter. If your transmitter is still not communicating with your pump, proceed to step 2. 2 Make sure your pump and your transmitter are located on the same side of your body to minimize any RF interference. If your transmitter is still not communicating with your pump, call your local representative for assistance.
Note: An alarm or alert occurs and a message appears when your transmitter has lost connection with your pump for 30 minutes.		

Storage and handling

Store the transmitter, charger, and tester in a clean, dry location at room temperature. If the transmitter is not in use, you must charge the transmitter at least once every 60 days.

CAUTION: Do not store the transmitter on the charger. If the transmitter is left on the charger for more than 60 days, the battery will be permanently damaged.

Disposal

Do not dispose the transmitter in unsorted municipal waste stream. Discard the transmitter according to local regulations for battery disposal, or contact your healthcare professional for disposal information.

Specifications

The essential performance (EP) of the transmitter is to measure and transmit to a monitoring device the sensing device's signal value(s) within the transmitter's accuracy requirements under the specified use conditions outlined in the system user guide and for the duration of the expected service life. If the transmitter experiences electromagnetic disturbances, either no or incorrect data may be transmitted. In such situations, refer to the operation, maintenance, and troubleshooting instructions within the applicable user guides. You may also use the tester to test if the transmitter is operating properly. If the transmitter is damaged or if it cannot communicate with the charger or pump, contact your local Medtronic support representative for assistance.

Biocompatibility	Transmitter: Complies with EN ISO 10993-1
Applied parts	Transmitter Sensor
Operating conditions	Transmitter temperature: 0°C to 45°C (32°F to 113°F) Caution: When operating the transmitter on a tester in air temperatures greater than 41°C (106°F), the temperature of the transmitter may exceed 43°C (109°F). Transmitter relative humidity: 10% to 95% with no condensation Transmitter pressure: 57.60 kPa to 106.17 kPa (8.4 psi to 15.4 psi) Charger temperature: 10°C to 40°C (50°F to 104°F) Charger relative humidity: 30% to 75% with no condensation
Storage conditions	Transmitter temperature: -20°C to 55°C (-4°F to 131°F) Transmitter relative humidity: up to 95% with no condensation Transmitter pressure: 57.6 kPa to 106 kPa (8.4 psi to 15.4 psi) Charger temperature: -10°C to 50°C (14°F to 122°F) Charger relative humidity: 10% to 95% with no condensation
Battery life	Transmitter: Seven days of continuous glucose monitoring immediately following a full charge. Charger: The charger uses one new AAA battery to charge the transmitter.
Transmitter frequency	2.4 GHz band, Bluetooth™* wireless technology (version 4.0)
Effective radiated power (ERP)	-12.05 dBm (0.06 mW)
Effective isotropic radiated power (EIRP)	-9.9 dBm (0.1 mW)
Operating range	Up to 1.8 meters (6 feet) in free-air
Transmitter expected service life	The transmitter expected service life is one year depending on patient usage.

Transmitter wireless communication

Quality of service

The transmitter and insulin pump connect via smart device connectivity. The transmitter sends glucose data and system-related alerts to the pump. The pump verifies the integrity of received data after wireless transmission.

Data security

The transmitter is designed to only accept radio frequency (RF) communications from recognized and linked devices. You must pair your pump with the transmitter before the pump will accept information from the transmitter.

The MiniMed™ insulin pumps and system components (meters and transmitters) ensure data security via proprietary means and data integrity using error checking processes, such as cyclic redundancy checks.

Traveling by air

Your transmitter is safe for use on commercial airlines. If questioned by airline personnel about the use of your device, please show them your Medical emergency card.

Guidance and manufacturer's declaration

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	CISPR 11 Group 1, Class B	The transmitter uses RF energy only for system communications. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. Note: The preceding statement is required by IEC 60601-1-2 for Group 1, Class B devices. Since the transmitter is battery powered, its emissions will not be affected by the establishment power supply and there is no evidence of any issues associated with the use of the system in domestic establishments.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
Immunity Test	IEC 60601-1-2:2014 Test Level	Max foreseeable use condition per IEC 60601-1-2:2014	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	For use in a typical domestic, commercial, or hospital environment.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
Immunity Test	IEC 60601-1-2:2014 Test Level	Max foreseeable use condition per IEC 60601-1-2:2014	Electromagnetic Environment Guidance
Conducted disturbances induced by RF fields	3 Vrms 150 kHz to 80 MHz 6 Vrms ISM bands between 150 kHz to 80 MHz	Not applicable	Requirement does not apply to this battery powered device.
Electrical fast transient/burst IEC 61000-4-4	±2 kV 100 kHz repetition frequency	Not applicable	Requirement does not apply to this battery powered device.
Surge IEC 61000-4-5	Line to Line: ±0.5 kV, ±1 kV Line to Ground: ±0.5 kV, ±1 kV, ±2 kV	Not applicable	Requirement does not apply to this battery powered device.
Note: U_T is the a.c. mains voltage prior to application of the test level.			
Voltage dips, short interruptions, and voltage variations on power supply lines IEC 61000-4-11	0% U_T ; 0.5 cycles (at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°) 0% U_T ; 1 cycle (at 0°) 70% for 25/30 cycles (at 0°) 0% for 250/300 cycles	Not applicable	Requirement does not apply to this battery powered device.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	For use in a typical domestic, commercial, or hospital environment.
Proximity fields from RF wireless communications equipment IEC 61000-4-3	IEC 60601-1-2:2014, Table 9	IEC 60601-1-2:2014, Table 9	For use in a typical domestic, commercial, or hospital environment.
Note: U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
Immunity Test	IEC 60601-1-2:2014 Test Level	Max foreseeable use condition per IEC 60601-1-2:2014	Electromagnetic Environment Guidance
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz	10 V/m 80 MHz to 6 GHz 80% AM at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of the transmitter than the recommended separation distance of 30 cm (12 in). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption, and reflection from structures, objects and people.			

Warranty

Medtronic MiniMed, Inc. (or such other legal entity as may be referred to as manufacturer on the labeling of this device "Medtronic MiniMed") warrants the Medtronic transmitter to the purchaser of the product against defects in material and workmanship for a period of one (1) year and the charger for up to one (1) year from the date of purchase.

During the warranty period, Medtronic MiniMed will replace or repair, at its discretion, any defective transmitter or charger, subject to the conditions and exclusions stated herein. This warranty applies only to new devices. In the event a transmitter or charger is replaced, the warranty period will not be extended past its original expiration date.

This warranty is valid only if the Medtronic transmitter or charger is used in accordance with the manufacturer's instructions. Without limitation, this warranty will not apply:

- If damage results from changes or modifications made to the transmitter or charger by the user, or third persons, after the date of purchase.

- If damage results from service or repairs performed by any person or entity other than the manufacturer.
- If damage results from a *Force Majeure* or other event beyond the control of the manufacturer.
- If damage results from negligence or improper use, including but not limited to: improper storage, submersion in water, physical abuse, (such as dropping).
- If damage results from use of the device in a manner other than according to the manufacturer's product labeling, instructions for use, or regulatory notifications.

This warranty shall be personal to the original purchaser. Any sale, rental or other transfer or use of the product covered by this warranty to or by a user other than the original purchaser shall cause this warranty to immediately terminate. This warranty does not apply to Glucose Sensors and other accessories.

The remedies provided for in this warranty are the exclusive remedies available for any breach hereof. Neither Medtronic MiniMed nor its suppliers or distributors shall be liable for any incidental, consequential, or special damage of any nature or kind caused by or arising out of a defect in the product.

All other conditions and warranties, other than mandatory statutory warranties, expressed or implied, are excluded, including the warranties of merchantability and fitness for a particular purpose.

This warranty gives the purchaser specific legal rights, and the purchaser may also have other rights that vary under local law. This warranty does not affect the purchaser's statutory rights.

Open Source Software (OSS) Disclosure


This document identifies the Open Source Software that may be separately called, executed, linked, affiliated, or otherwise utilized by this product.

Such Open Source Software is licensed to users subject to the terms and conditions of the separate software license agreement for such Open Source Software.

Use of the Open Source Software by you shall be governed entirely by the terms and conditions of such license.

The source/object code and applicable license for the Open Source Software can be obtained at the following site: <http://www.ouah.org/ogay/hmac/>.

Icon Table

	Serial number
	Catalogue or model number
(1x)	One transmitter, charger, and serter per container/package
(2x)	Two testers per container/package
	Date of manufacture (YYYY-MM-DD)
	Manufacturer
	Must refer to instruction manual before every use (appears blue on label).
	Temperature limit
	Non-ionizing electromagnetic radiation
	Configuration or unique version identifier
	Degree of protection against electric shock: Type BF applied part
IP48	Transmitter: 4 is the level of protection against solid objects with a diameter above 1mm. 8 is the level of protection against the effects of continuous immersion in water [2.4 meters (8 feet) immersion for 30 minutes].
	Humidity limitation
CE 0459	Marking of Conformity: This symbol means the device fully complies with MDD 93/42/EEC (NB 0459).
	Fragile, handle with care
	Keep dry
	Recycle cardboard, paper, plastic packaging supplies and unwanted written material.

	WEEE Initiative: DO NOT THROW IN TRASH. Recycle device according to local disposal requirements.
	Magnetic Resonance (MR) Unsafe: keep away from magnetic resonance imaging (MRI) equipment
	Recharge by date
	Bluetooth™ wireless technology
FCC ID: OH27821	Complies with United States Radio communication requirements.
IC: 3408B-7821	Complies with Innovation, Science and Economic Development Canada Radio communication requirements.
	Authorized representative in the European community
	Complies with Australian Radio communications requirements.
R _x Only	Prescription only
	Marking of Conformity: This symbol means the device fully complies with MDD 93/42/EEC.
	Consult instructions for use
	Caution: consult instructions for use for important warnings or precautions not found on the label

© 2020 Medtronic MiniMed, Inc. All rights reserved.

Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. ™* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

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ATTESTATION STATEMENT

As a responsible member of Medtronic MiniMed, I hereby attest that I believe the Appendix to the IFU attached herein, to be a true and accurate copy.

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
On July 16th 2021 before me, Kusum Manuwani Notary Public
(date) (insert name and title of the officer)

personally appeared Katherine Rojas, who proved to me on the basis of satisfactory evidence to be the person(~~s~~) whose name(~~s~~) is/~~are~~ subscribed to the within instrument and acknowledged to me that ~~he~~/she/~~they~~ executed the same in ~~his~~/her/~~their~~ authorized capacity(~~ies~~), and that by ~~his~~/her/~~their~~ signature(~~s~~) on the instrument the person(~~s~~), or the entity upon behalf of which the person(~~s~~) acted, executed the instrument. I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

WITNESS my hand and official seal.

Signature Kusum Manuwani (Seal)

Guardian Link Transmitter Kit (3) for MiniMed Pumps with Continuous Glucose Monitoring System

Appendix to IFU

Name of medical device	Guardian Link Transmitter Kit (3) for MiniMed Pumps with Continuous Glucose Monitoring System. Guardian Link Transmitter Kit (3) for MiniMed Pumps with Continuous Glucose Monitoring System: I. Set I (MMT-7810W3), consisting of: 1. Guardian Link Transmitter (3) (MMT-7811WW) – 1 p.; 2. Charger (MMT-7715) – 1 p. (if necessary) 3. Universal tester (MMT-7736L) – 2 p.; (if necessary) 4. One-press Sertter (MMT-7512) (if necessary) 5. Operational documentation: II. Set II (MMT-7811W3), consisting of: 1. Guardian Link Transmitter (3) (MMT-7811WW) – 1 p.; 2. Operational documentation: III. Guardian Link Transmitter Kit (3) with Bluetooth for MiniMed Pumps with Continuous Glucose Monitoring System: I. Set I (MMT-7910W3), consisting of: 1. Guardian Link Transmitter (3) (MMT-7911WW) – 1 p.; 2. Charger (MMT-7715) – 1 p.; (if necessary) 3. Universal tester (MMT-7736L) – 2 p.; (if necessary) 4. One-press Sertter (MMT-7512) (if necessary) 5. Operational documentation: IV. Set II (MMT-7911W3), consisting of: 1. Guardian Link Transmitter (3) (MMT-7911WW) – 1 p.; 2. Operational documentation Glossary <table border="1"> <thead> <tr> <th>Key term</th><th>Additional terms that may be found in the operating instructions</th></tr> </thead> <tbody> <tr> <td>Guardian Link Transmitter (3)</td><td>The Guardian™ Link (3) transmitter with Bluetooth™ wireless technology, device, transmitter</td></tr> <tr> <td>Continuous glucose monitoring (HMT)</td><td>Continuous glucose monitoring (CGM)</td></tr> <tr> <td>One-press Sertter</td><td>One-press sertter</td></tr> </tbody> </table>	Key term	Additional terms that may be found in the operating instructions	Guardian Link Transmitter (3)	The Guardian™ Link (3) transmitter with Bluetooth™ wireless technology, device, transmitter	Continuous glucose monitoring (HMT)	Continuous glucose monitoring (CGM)	One-press Sertter	One-press sertter
Key term	Additional terms that may be found in the operating instructions								
Guardian Link Transmitter (3)	The Guardian™ Link (3) transmitter with Bluetooth™ wireless technology, device, transmitter								
Continuous glucose monitoring (HMT)	Continuous glucose monitoring (CGM)								
One-press Sertter	One-press sertter								
Prescription of medical device	Intended for use by single-patient or multiple-patient use for processing and storing glucose sensor values, as well as for transmitting this data to compatible Medtronic MiniMed devices.								
Indications and Contraindications for medical device	<u>The Guardian Link (3) transmitter</u> Indications for use								

• Available on request

The Guardian Link (3) transmitter is intended for use with the Guardian Sensor (3). The transmitter powers the glucose sensor, collects and calculates sensor data, and wirelessly sends the data to a pump. The transmitter is intended for single-patient or multiple-patient use.

Contraindications

None known.

The Guardian Link (3) Transmitter (MMT-7911)**Indications for use**

The Guardian Link (3) transmitter (MMT-7911) is a rechargeable device and powers the glucose sensor, collects and calculates the sensor data, and sends the data to a compatible MiniMed insulin pump system with smart device connectivity for the management of diabetes mellitus. The transmitter is only compatible with the Guardian Sensor (3) glucose sensor (MMT-7020) and is indicated for single-patient or multiple-patient use.

Contraindications

None known.

Serter**Indications for use**

The serter is used as an aid for inserting the sensor. It is indicated for single-patient use and is not intended for multiple patient use.

Contraindications

Do not use the serter on products other than the approved Medtronic MiniMed glucose sensors. Medtronic cannot guarantee the safety or efficacy of this product if used with other products.

Tester**Indications for use**

The Watertight Tester is an accessory item included with select Medtronic Diabetes continuous glucose transmitters. The tester connects to a transmitter and is used to test the transmitter. The tester is also used to form a waterproof seal to protect the connector on the transmitter while the transmitter is being cleaned.


Contraindications

None known.

Charger**Indications for use**

The charger is an accessory item included with select Medtronic continuous glucose recorder and transmitter devices. It is used to recharge the internal battery for these devices, as needed.


* Available on request

	Contraindications None known.	
Information on potential consumers of medical device, use conditions	Intended for use by patient or physician. Intended for use by adult and children both in inpatient and home conditions.	
Description of the main functional parts of medical device	The tolerances are $\pm 10\%$	
	Transmitter specifications (MMT-7811)	
		
	Protection class	IP48
	Type of applied part	BF
	Rating	No ratings (Transmitter Internally Powered 3.7Vdc to be charged only via Charger) Lithium Battery: Rated 3.7 V, 25 mAH (Max maximum charging current 75mA)
	Electrical shock protection class	Class II
	Software Version	V2.1A
	Software safety class	C
	Software release date	August 25, 2015
	Dimensions, W × L × H	3.5 cm (+0.018/-0.008 cm) (w) × 4.1 cm (+0.018/-0.008 cm) (l) × 1.3 cm (± 0.013 cm) (h)
	Mass	< 14.2 grams
	Sterilization	Non-sterile device
	Shelf life	Six (6) months
	Operational life	One (1) year

• Available on request

Battery life	Seven (7) days following a full charge
Battery	Lithium polymer
Transmitter frequency	2.4 GHz, 2M60G1D modulation, less than 1 mW ERP
Maximum output power (EIRP)	1 mW (0 dBm)
Radio Frequency (RF) communications	Pump to transmitter frequency: 2.4 GHz; proprietary Medtronic protocol; range up to 1.8 meters (6 feet) Utilizes the IEEE 802.15.4 protocol with proprietary data format
Operating mode	Continuous
The maximum permissible time for establishing the operating mode, calculated from the moment of switching on or starting the product	There is no duty cycle for the operating mode as it is continuous.
Sterilization method (s)	Not sterile, designed for routine cleaning
Suitability for use in high oxygen environments	Not intended for use in high oxygen environments.

Transmitter specifications (MMT-7911)



Protection class	IP48
Type of applied part	BF
Ratings	No ratings (Transmitter Internally Powered 3.7Vdc to be charged only via Charger) Lithium Battery: Rated 3.7 V, 25 mAH (Max maximum charging current 75mA)
Electrical shock protection class	Class II

• Available on request

Software Version	V1.0A
Software safety class	C
Software release date	October 18, 2019
Dimensions, W × L × H	2.84 cm (w) × 3.57 cm (l) × 0.93 cm (h)
Mass	4.75 grams
Sterilization	Non-sterile device
Shelf life	Six (6) months
Operational life	One (1) year
Battery life	Seven (7) days following a full charge
Battery	Lithium polymer
Transmitter frequency	2.4 GHz, 2M60G1D modulation, less than 1 mW ERP
Maximum output power (EIRP)	-9.9 dBm (0.1 mW)
Radio Frequency (RF) communications	Bluetooth Low Energy (BLE, 2.4 GHz) wireless communication protocol
Operating mode	Continuous
The maximum permissible time for establishing the operating mode, calculated from the moment of switching on or starting the product	There is no duty cycle for the operating mode as it is continuous.
Sterilization method (s)	Not sterile, designed for routine cleaning
Suitability for use in high oxygen environments	Not intended for use in high oxygen environments.
Tester specifications	


• Available on request

Dimensions, L × W × H	2.98 cm (± 0.013 cm) × 1.94 cm (± 0.013 cm) × 0.7 cm (± 0.013 cm) [1.5 inches × 0.8 inches × 0.3 inches]
Mass, maximum	3 grams (0.11 oz)

Charger specifications

* Available on request


Dimensions, H × W × L	2.9 cm (± 0.013 cm) × 4.0 cm (± 0.013) × 6.1 cm (+0.015/-0.01) (1.75 x 2.5 x 1.15 inches)
Mass, maximum	25 g without AAA battery
Battery required	AAA alkaline battery
Battery capacity	1250 mAh
Nominal battery voltage	1,5 V
Protection class	IPX0

• Available on request

One-press Serter specifications

Available on request


Dimensions, H × W × L	7.85 cm (± 0.013 cm) x 6.91 cm (± 0.025 cm) x 5.77 cm (± 0.025 cm) (3.25 in x 2.50 in x 3.00 in)
Mass	112.5 g +/- 15%
Operating conditions	The Serter shall operate within a normal room temperature of 10 °C to 40 °C
Storage conditions	Temperature: -20 °C +0/-2 °C to 55 °C +2/-0 °C. Atmospheric pressure: 49.64 kPa +0/-3.5 kPa to 106.18 kPa +3.5/-0 kPa (7.2 psia +0/-0.5 psia to 15.4 psia +0.5/-0 psia). Relative humidity: 15 % to 85 % condensing
Shelf life	There is no shelf life requirement associated with the Serter
Operating life	The serter has a minimum service life of one year

Description of all accessories, medical devices or non-medical devices, intended for use with the medical device submitted for registration

The following is a brief description of the units and compatible devices to the Guardian Link (3) Transmitter:

Charger (MMT-7715)

Available on request

1 - PCBA
2 - Housing, Top
3 - Housing, Bottom
4 - Battery Door
5 - Connector
6 - Label, Top
7 - Label, Bottom
9 - Contact, Positive
10 - Contact, Negative
11 - Connector


The charger (MMT-7715) is an accessory device intended to charge the rechargeable lithium polymer battery contained within the Guardian Link (3) Transmitter (MMT-7811). The Charger has been classified in accordance to Annex IX of MDD (93/42/EEC) as a Class I device, under Rule 12.

• Available on request

The transmitter is charged automatically following connection to the charger. The charger automatically terminates charging upon full charge.

The charger is powered by an internal, replaceable, commercially available AAA alkaline battery and is designed with sliding rear cover to facilitate changing the battery. The charger uses visual indicators (green/red LEDs) to show the status of transmitter device charging. The indicator light also communicates any problem with the transmitter device battery, or the charger battery.

Universal tester (MMT-7736L)

3 – O-Ring
4 – O-Ring
5 – Dual Trace Circuit
8 – Test Plug

Universal tester (MMT-7736L) is an accessory device intended to facilitate troubleshooting and the cleaning of Guardian Link (3) Transmitter (MMT-7811).

The Tester is considered a non-active medical device and has been classified in accordance with Annex IX of MDD (93/42/EEC) as a Class I device, under Rule 1.

The Tester (MMT-7736L) is assembled from a plastic housing and a printed circuit board (PCB). The tester is designed with a molded handle and plug that mate with transmitter device connector in the same manner that a sensor would be attached to the transmitter device.

Available on request

The tester operates as a sensor simulator, creating a signal current at a level within a range of an *in vivo* sensor during normal operation. The tester also provides a watertight seal for cleaning the transmitter. The tester is used as a handle to aid in holding the transmitter during cleaning or troubleshooting. The tester PCB contains a resistor. Once the tester is plugged into the transmitter device, the tester connects the device counter and reference connections to cause the transmitter device to start operating.

Compatible Devices

The Guardian Link (3) Transmitter is intended to be used as a component of select Medtronic MiniMed continuous glucose monitoring or sensor enabled pump systems, and is compatible with following devices:

Model Number	Device Name
MMT-7020	Guardian (3) Sensor
MMT-1780	MiniMed 670G Insulin Pump
MMT-1711, MMT-1712, MMT-1751, MMT-1752	MiniMed 640G Insulin Pump

Receiving device: MiniMed 670G Insulin Pump, MiniMed 640G Insulin Pump

The MiniMed 670G Insulin Pump, MiniMed 640G Insulin Pump are an ambulatory, battery operated, rate-programmable micro-infusion pump. In addition, to insulin delivery, the MiniMed Pump, is designed for use with Medtronic glucose sensing systems. It is designed to receive, display, and store real-time sensor glucose values received from a compatible transmitter. The sensor signals from the sensor are converted by the compatible transmitter to sensor glucose values (based on calibration values from a blood glucose meter) which are transmitted to the pump for display to the user. The glucose sensor values are transmitted via radiofrequency (RF) telemetry every 5 minutes.

The transmitter is capable of transmitting and receiving communication from compatible Medtronic devices (i.e., Medtronic MiniMed 600G series pump) via a Medtronic proprietary communication protocol within the 2.4 – 2.4835 GHz ISM band. The chip controlling the 2.4 GHz radio contains a modified version of the IEEE 802.15.4 protocol stack. The protocol enables the formation of a network and two-way communication between devices.

The transmitter and all compatible system devices utilize the same Medtronic proprietary telemetry and are linked in a star IEEE 802.15.4 based beacon-enabled network configuration, in which the pump serves as the coordinator and the other devices in the network as end-nodes. The coordinator (i.e. 600G Series Insulin Pump) has the ability to determine which channel has the least amount of RF noise, and if necessary, it will change to that channel. The end-nodes will follow the coordinator and the network will be re-established in this new channel.

Under the beacon-enabled network scheme, the network coordinator broadcasts a periodic beacon then turns the radio receiver ON for about 15 ms to allow end-nodes to transmit any packets to the coordinator. In order to assure network synchronization, each

• Available on request


	end-node turns on its receiver and listens for the coordinator beacon on a specific programmable interval. When communications are initiated by the pump, a "request" to communicate with one or more end-nodes is posted on the beacon. When communications are initiated by an end-node, the end-node sends the information immediately after a beacon. The System network is a Wireless Body Area Network (WBAN) and any single user with a compatible MiniMed 600G series pump can be concurrently connected with up to seven end-nodes. All devices in the network are uniquely identified by their Media Access Control Identifier (MACID). Using a commercially available radio, the TRANSMITTER communicates in the 2.4 GHz ISM band (2.4 – 2.4835 GHz) using a proprietary protocol based on IEEE 802.15.4. Radio telemetry (Tel-D) requirements. Compatible Sensor: Guardian (3) Sensor The sensor continuously converts tiny amounts of glucose from the interstitial fluid under the skin, into an electronic signal. The CGM system then uses these signals to provide sensor glucose values. The glucose sensing principle of the glucose sensor is based on the oxidation of glucose. The coating that is placed on the electrode surface contains the enzyme, glucose oxidase. This enzyme catalyzes the reaction between glucose and oxygen to produce an electrical signal
Risks of using a medical device, expected and predictable side effects associated with the intended use of a medical device	Common risks associated with using the transmitter include: There may be minor skin irritation when the product is in contact with skin.
Methods and apparatus for disinfection and sterilization cleaning	Cleaning the transmitter The transmitter is intended for personal use at home (single-patient use) or for use in healthcare facilities (multiple-patient use). Single-patient use requires cleaning after each use, while multiple-patient use requires cleaning and disinfection after each use. When using the transmitter in a healthcare facility, always follow the cleaning and disinfecting procedure for multiple-patient use. Warning: Do not discard the transmitter in a medical waste container or expose it to extreme heat. The transmitter contains a battery that may ignite, and result in serious injury. <i>Note: The tester is a required component for cleaning the transmitter.</i> For single-patient use Always clean the transmitter after each use. To clean the transmitter, you need the following materials: • mild liquid soap (the device was validated with dish washing liquid, but any commercially available mild liquid soap is acceptable) • soft-bristled toddler toothbrush • container • clean, lint-free dry cloths Use life The transmitter can be cleaned up to 122 times or one year, whichever comes first. Discard the transmitter at this point. If you continue to use the transmitter beyond 122 times or one year, the cleaning process may damage the device. Contact Medtronic to order a new transmitter.

• Available on request

Warning: Do not use the device if you see any cracking, flaking, or damage to the housing. Cracking, flaking, or damage to the housing are signs of deterioration. Deterioration of the housing can affect the ability to properly clean the transmitter, and result in serious injury. Call your local representative and discard the device according to local regulations for battery disposal (non-incineration), or contact your healthcare professional for disposal information.


To clean the transmitter:

- 1 Wash your hands thoroughly.
- 2 Attach the tester to the transmitter.

- 3 If there is adhesive residue on the transmitter, see Removing adhesive residue in IFU.

- 4 Rinse the transmitter under room temperature tap water for at least one minute, and until visibly clean. Make sure all hard-to-reach areas are rinsed completely.

- 5 Prepare a mild liquid soap solution using 5 milliliters (1 teaspoon) of mild liquid soap per 3.8 liters (1 gallon) of room temperature tap water.

- 6 With the tester still attached, submerge the transmitter in the mild liquid soap solution and soak for one minute.


• Available on request

7 Holding the tester, brush the entire surface of the transmitter using a soft-bristled toddler toothbrush. Make sure to brush all hard-to-reach areas until visibly clean.

8 Rinse the transmitter under running room temperature tap water for at least one minute, and until all visible liquid soap is gone.

9 Dry the transmitter and tester with a clean, dry cloth.

10 Place the transmitter and tester on a clean, dry cloth and air dry them completely.

11 Disconnect the tester from the transmitter by gently squeezing the arms of the tester.

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Warning: You must adhere to Standard Precautions when handling or using this device. All parts of the system should be considered potentially infectious and are capable of transmitting blood-borne pathogens between patients and healthcare professionals.

The transmitter must be disinfected after use on each patient. This device may only be used for testing multiple patients when Standard Precautions and Medtronic's disinfection procedures are followed.

To clean and disinfect the transmitter, you need the following materials:

- gloves
- mild liquid soap (the device was validated with dish washing liquid, but any commercially available mild liquid soap is acceptable)
- soft-bristled toddler toothbrush
- 8.25% bleach
- two containers
- clean, lint-free dry cloths

Use life

The transmitter can be cleaned up to 122 times or one year, whichever comes first.

Discard the transmitter at this point. If you continue to use the transmitter beyond 122 times or one year, the cleaning process may damage the device. Contact Medtronic to order a new transmitter.

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
- 1 Wash your hands and put on gloves.
- 2 Inspect the inside of the transmitter connector opening for any sign of body fluid.


For instructions on how to inspect the connector pins, see Inspecting the transmitter connector pins, on page 8.

Caution: The person inspecting the transmitter must have sufficient vision that enables him or her to see small drops of body fluid or debris.

Warning: If you see any body fluid in the connector opening, you must discard the transmitter. Because the transmitter contains a battery, do not discard in a bio-waste container. Instead, continue to clean and disinfect the transmitter, and then discard according to local regulations for battery disposal (non-incineration).

- 3 Attach the tester to the transmitter.

4 If there is adhesive residue on the transmitter, see Removing adhesive residue in IFU.

5 Rinse the transmitter under room temperature tap water for at least one minute, and until visibly clean. Make sure all hard-to-reach areas are rinsed completely.

6 Prepare a mild liquid soap solution using 5 milliliters (1 teaspoon) of mild liquid soap per 3.8 liters (1 gallon) of room temperature tap water. Make sure to prepare a fresh solution for each use.

7 With the tester still attached, fully submerge the transmitter in the mild liquid soap solution and soak for one minute.

8 Holding the tester, brush the entire surface of the transmitter using a soft-bristled toddler toothbrush. Make sure to brush all hard-to-reach areas until visibly clean.

9 Rinse the transmitter under running room temperature tap water for at least one minute, and until all visible liquid soap is gone.

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11 Prepare a 1:10 bleach solution by using one (1) part 8.25% bleach to nine (9) parts water, for a final concentration of 0.8%. Make sure to prepare a fresh solution for each use.

12 Make sure that you have completed the previous cleaning steps before disinfection. With the tester still attached, soak the transmitter in the bleach solution for 20 minutes.

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
Warning If you saw any body fluid inside the connector opening on earlier inspection, you must now discard the transmitter with tester still attached, according to local regulations for battery disposal (non incineration).

15 Disconnect the tester from the transmitter by gently squeezing the arms of the tester.

16 Inspect the housing of the transmitter for any signs of cracking, flaking, or damage. If you see any of these signs, you must now discard the disinfected transmitter according to local regulations for battery disposal (non-incineration).

Warning: Do not use the device if you see any cracking, flaking, or damage to the housing. Cracking, flaking, or damage to the housing are signs of deterioration. Deterioration of the housing can affect the ability to properly clean the transmitter, and result in serious injury. Call your local representative and discard the device according to local regulations for battery disposal (non-incineration), or contact your healthcare professional for disposal information.

17 Discard the used gloves and thoroughly wash hands with soap and water.

Removing adhesive residue

You may need to perform this procedure if there is adhesive residue present on the transmitter. If you visually inspect the transmitter and see adhesive residue on it, follow the instructions below.

To remove adhesive residue, you need medical adhesive remover (such as Detachol® - which is a mineral spirit) and cotton swabs.

Note: During testing, Medtronic MiniMed used Detachol to remove the adhesive residue from the transmitter. Detachol is recommended for use but may not be available in all countries.

To remove adhesive residue:

1 Make sure the tester is attached to the transmitter.

2 Holding the tester, saturate a cotton swab in the Detachol solution and gently rub the adhesive residue on the transmitter until it is fully removed.

3 Continue with the appropriate single-patient or multiple-patient cleaning procedure.

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Caution: Do not immerse the charger in water or any other cleaning agent. The charger is not waterproof. Water can damage the charger, and cause the device to malfunction.

Warning: Dispose the charger according to the local regulations for battery disposal, or contact your healthcare professional for disposal information. The charger may ignite upon incineration.

Warning: For multiple-patient use, always clean and disinfect the transmitter after removing it from the patient and before attaching it to the charger. If blood comes in contact with any surface areas of the charger, the contaminated device must be discarded. The charger contains a battery which may ignite upon incineration.

To clean the charger:

1 Wash your hands thoroughly.

2 Use a damp cloth with mild cleaning solution, such as a dishwashing detergent, to clean any dirt or foreign material from the outside of the charger. Never use organic solvents, such as paint thinner or acetone, to clean the charger.

3 Place the charger on a clean, dry cloth and air dry for two to three minutes

Utilization

Discard the transmitter according to local regulations for battery disposal, or contact your healthcare professional for disposal information.

Medical waste class for unused device: A

Medical waste class for used device: B

Medical waste class of electronic components: F

Warranty

Manufacturer guarantees compliance of the device with standards and product requirements subject to instructions of use, transport, storage and utilization regulations.

The list of national and international regulatory documents/standards, which matches medical device

Standard	Description
BS EN ISO 13485:2016	Medical devices – Quality management systems – Requirements for regulatory purposes
BS EN ISO 14971:2012	Medical devices – Application of risk management to medical devices (ISO 14971:2007, Corrected version 2007-10-01)

Available on request

BS EN 62304:2006+A1: 2015	Medical device software – Software life-cycle processes
ETSI EN 300 328 v2.1.1: 2016	Electromagnetic compatibility and Radio spectrum Matters (ERM); Wideband transmission systems; Data transmission equipment operating in the 2,4 GHz ISM band and using wide band modulation techniques; Harmonized EN covering the essential requirements of article 3.2 of the R&TTE Directive
ETSI EN 301 489-1 v2.1.1:2017	Electromagnetic compatibility and Radio spectrum Matters (ERM); Electro Magnetic Compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements
ETSI EN 301 489-17 v3.1.1: 2017	Electromagnetic compatibility and Radio spectrum Matters (ERM); Electro Magnetic Compatibility (EMC) standard for radio equipment; Part 17: Specific conditions for Broadband Data Transmission Systems
BS EN 1041:2008+A1: 2013	Information supplied by the manufacturer of medical devices - AMD: October 31, 2013
BS EN ISO 15223-1:2016	Medical devices – Symbols to be used with medical device labels, labelling, and information to be supplied – Part 1: General requirements.
ASTM D4169-16:2016	Standard Practice for Performance Testing of Shipping containers and system
BS EN ISO 10993-1:2009/C:2010	Biological evaluation of medical devices Part 1: Evaluation and testing within a risk management process - CORR: July 31, 2010; CORR: November 30, 2010
ISO 10993-1:2018	Biological evaluation of medical devices Part 1: Evaluation and testing within a risk management process
BS EN ISO 10993-5:2009	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
BS EN ISO 10993-10: 2013	Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
BS EN ISO 10993-12: 2012	Biological evaluation of medical devices- Part 12: Sample preparation and reference materials
BS EN ISO 10993-17: 2009	Biological evaluation of medical devices- Part 17: Establishment of allowable limits for leachable substances
BS EN ISO 10993-18: 2009	Biological evaluation of medical devices- Part 18: Chemical characterization of materials
BS EN ISO 14155:2011	Clinical investigation of medical devices for human subjects - Good clinical practice - CORR: September 30, 2011
BS EN 60601-1:2006/+A1: 2013	Medical electrical equipment – Part 1: General requirements for basic safety and essential performance (IEC 60601-1:2005/A1:2012)
BS EN 60601-1-2: 2015	Medical electrical equipment – Part 1-2: General requirements for safety – Collateral standard: Electromagnetic compatibility – Requirements and tests
BS EN 60601-1-6: 2010	Medical electrical equipment – Part 1-6: General requirements for safety – Collateral standard: Usability
BS EN 60601-1-11: 2010	Medical electrical equipment - Part 1-11: general requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment

• Available on request

	BS EN 62479: 2010	Assessment of the compliance of low power electronic and electrical equipment with the basic restrictions related to human exposure to electromagnetic fields (10 MHz to 300 GHz)
	IEC 60529:2013	Degrees of Protection Provided by Enclosures (IP Code)
	BS EN 62366: 2008	Medical devices – Application of usability engineering to medical devices
	BS EN 62366-1:2015	Medical devices – Part 1: Application of usability engineering to medical devices
	BSI BS EN 17664: 2017	Processing of health care products. Information to be provided by the medical device manufacturer for the processing of medical devices. Part 1: Devices that require cleaning followed by disinfection and/or sterilization
	AAMI TIR12: 2010	Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers
	AAMI TIR30: 2011	A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices – FDA RECOGNIZED
	Japanese Radio Law (JRL) Notice 88	Ordinance Concerning Technical Regulations Conformity Certification of Specified Radio Equipment
	RSS-247:2015	Low-power License-exempt Radio communication Devices (All Frequency Bands): Category I Equipment
	47 CFR Part 15	RADIO FREQUENCY DEVICES
	RTCA DO-160G	Environmental Conditions and Test Procedures for Airborne Equipment
	AAMI TIR69: 2017	Risk Management of Radio-Frequency Wireless Coexistence for Medical Devices and Systems
	<i>MMT-7911 Only</i>	
	Bluetooth Specification v4.0	Bluetooth Core Specification Version 4.0
Customer Service	Contact for all questions, including complaints, repair and service of medical devices the Authorized Representative in Russia: Medtronic, LLC, 123112, Russia, Moscow, Presnenskaya nab, h. 10, fl.9, pom. III, ap.41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru	

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