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CE 0192

1. Indications for Use

DIA-ROOT BIO Sealer is a MTA (mineral trioxide aggregate) based root canal sealer that provides complete and permanent sealing of root canals. It can be used with or without root canal obturation materials.

2. Instructions

1. Mix the product with water in a 1:1 ratio.
2. Apply the mixture to the root canal using a syringe or a paper point.
3. Allow the mixture to set for 10 minutes.
4. Remove the excess material with a paper point.
5. The mixture is ready for use.

3) If the contents of the product come into contact with skin, immediately flush with plenty of water.

4) Do not use the product after the expiration date.

5) Do not store or use the contents in any container except the one provided by the manufacturer.


C. Contraindications

1) Do not use for purpose other than filling of root canals for dental purposes.

2) In case of allergic reaction, discontinue use immediately and consult a physician.

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horns with a ne

Mix DIA-ROOT BIO MTA with distilled water
 at the following ratio and add distilled water
 until the mixture hardens or dries according to the
 temperature and humidity. Fill the cavity
 and mix.
 * Mixing Ratio = 0.5g : 0.225g (0.225 cc
 (DIA-ROOT BIO MTA : Distilled Water)
 Remove excess water with sterilized cotton
 and compactly squeeze the empty space inside.
 Before DIA-ROOT BIO MTA hardens, add
 DIA-ROOT BIO MTA to fill the cavity.
 After setting is complete, apply a strong
 pressure to the hardened DIA-ROOT BIO MTA
 and check for...
 ...entirely and proceed
 ...as retreatment may cause pulpal
 ...the top part of DIA-ROOT BIO MTA
 ...apply resin or glass ionomer according to...
 ...partial pulpotomy of posterior teeth)
 1) Isolate the tooth with rubber dam and form
 a high speed bur and avoiding tissues close
 2) Remove decay closest to pulp thoroughly
 high speed diamond bur.
 3) Open and irrigate all pulp horns with a ne...

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