HYCARE
INTERNATIONAL CO., LTD

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E-Mail : info@hycare-int.com

Handwritten signature

INSTRUCTION FOR USE FOR MEDICAL DEVICE:
Medical diagnostic (examination) disposable non-sterile
gloves "BENOVY" natural latex

Managing director

HYCARE
INTERNATIONAL CO., LTD

(signature)

/Mr. Kevin Gostner/

2019

No. 161007 25 AUG 2021

BKK100

Certified Genuine Signature(s) of

KEVIN GOSTNER

(MR. EK-ATT THITARAM)

Minister-Counsellor

Ministry of Foreign Affairs of Thailand

(The Ministry of Foreign Affairs assumes no
responsibility for the contents of the documents)

3801123

TABLE OF CONTENTS

1 Data on the medical device	3
1.1 Name of the medical device	3
1.2 Manufacturer (producer) of the medical device.....	5
1.3 Purpose of the medical device	5
1.4 Scope of the medical device.....	5
1.5 Intended users of the medical device	5
1.6 Application conditions of the medical device.....	5
1.7 Type of contact with human body.....	5
1.8 Principle of operation of the medical device	6
1.9 Indications for use of the medical device.....	7
1.10 Contradictions of the medical device.....	7
1.11 Shelf life of the medical device.....	7
2 Warning and precautions	8
3 Instruction for use of the medical device	9
4 Transportation, storage and disposal conditions of the medical device.....	10
4.1 Transportation conditions of the medical device	10
4.2 Storage conditions of the medical device	10
4.3 Disposal conditions of the medical device.....	10
5 Cleaning and disinfection methods of the medical device.....	10
6 Specification of the medical device	11
6.1 Description of the medical device.....	11
6.2 Physical specification of the medical device	13
6.2.1 Physical and mechanical performance of the medical device.....	13
6.2.2 Inspection levels and acceptable quality level of the medical device.....	13
7 Adverse reactions.....	14
8 Circumstances, in which the consumer must consult healthcare professional	14
9 Manufacturer's warranties	14
10 Data on Authorized representative of the manufacturer	14
11 Scope of delivery of the medical device	14
12 List of national standards, applied by manufacturer (producer) of the medical device	15

TABLE OF CONTENTS

1 Data on the medical device	3
1.1 Name of the medical device	3
1.2 Manufacturer (producer) of the medical device	5
1.3 Purpose of the medical device	5
1.4 Scope of the medical device	5
1.5 Intended users of the medical device	5
1.6 Application conditions of the medical device	5
1.7 Type of contact with human body	5
1.8 Principle of operation of the medical device	6
1.9 Indications for use of the medical device	7
1.10 Contradictions of the medical device	7
1.11 Shelf life of the medical device	7
2 Warning and precautions	8
3 Instruction for use of the medical device	9
4 Transportation, storage and disposal conditions of the medical device	10
4.1 Transportation conditions of the medical device	10
4.2 Storage conditions of the medical device	10
4.3 Disposal conditions of the medical device	10
5 Cleaning and disinfection methods of the medical device	10
6 Specification of the medical device	11
6.1 Description of the medical device	11
6.2 Physical specification of the medical device	13
6.2.1 Physical and mechanical performance of the medical device	13
6.2.2 Inspection levels and acceptable quality level of the medical device	13
7 Adverse reactions	14
8 Circumstances, in which the consumer must consult healthcare professional	14
9 Manufacturer's warranties	14
10 Data on Authorized representative of the manufacturer	14
11 Scope of delivery of the medical device	14
12 List of national standards, applied by manufacturer (producer) of the medical device	15

1 Data on the medical device

1.1 Name of the medical device

«Medical diagnostic (examination) disposable non-sterile gloves “BENOVY” natural latex», in hardware versions:

1. Medical diagnostic (examination) disposable non-sterile gloves “BENOVY” natural latex smooth powdered (beaded cuff), color: natural, blue, pink, yellow, size: XS, S, M, L, XL. Series: 5-P. Type 122560.
2. Medical diagnostic (examination) disposable non-sterile gloves “BENOVY” natural latex textured powder free with inner polymer coating (beaded cuff), color: natural, size: XS, S, M, L, XL. Series: 55-1P. Type 122540.
3. Medical diagnostic (examination) disposable non-sterile gloves “BENOVY” natural latex textured powder free with inner and outer polymer coating (Low Risk Allergy) (beaded cuff), color: orange, green mint, light pink, lilac, fuschia, blue, black, light green, yellow, size: XS, S, M, L, XL. Series: 55-2P. Type 122540.
4. Medical diagnostic (examination) disposable non-sterile gloves “BENOVY” natural latex textured powder free with inner polymer coating with Aloe Vera extract (beaded cuff), color: natural, size: XS, S, M, L, XL. Series: 55-1P/A. Type 122540.
5. Medical diagnostic (examination) disposable non-sterile gloves “BENOVY” natural latex textured powder free with inner and outer polymer coating (Low Risk Allergy) with Aloe Vera extract (beaded cuff), color: orange, green mint, light pink, lilac, fuschia, blue, black, light green, yellow, size: XS, S, M, L, XL. Series: 55-2P/A. Type 122540.
6. Medical diagnostic (examination) disposable non-sterile gloves “BENOVY” natural latex textured powder free with inner polymer coating (beaded cuff), color: natural, size: XS, S, M, L, XL. Series: 62-1P. Type 122540.
7. Medical diagnostic (examination) disposable non-sterile gloves “BENOVY” natural latex textured powder free with inner and outer polymer coating (Low Risk Allergy) (beaded cuff), color: orange, green mint, light pink, lilac, fuschia, cornflower blue, black, light green, yellow, size: XS, S, M, L, XL. Series: 62-2P. Type 122540.
8. Medical diagnostic (examination) disposable non-sterile gloves “BENOVY” natural latex textured powder free with inner polymer coating with Aloe Vera extract (beaded cuff), color: natural, size: XS, S, M, L, XL. Series: 62-1P/A. Type 122540.
9. Medical diagnostic (examination) disposable non-sterile gloves “BENOVY” natural latex textured powder free with inner and outer polymer coating (Low Risk Allergy) with Aloe Vera extract (beaded cuff), color: orange, green mint, light pink, lilac, fuschia, blue, black, light green, yellow, size: XS, S, M, L, XL. Series: 62-2P/A. Type 122540.
10. Medical diagnostic (examination) disposable non-sterile gloves “BENOVY” natural latex textured powder free chlorinated (beaded cuff), color: natural, size: XS, S, M, L, XL. Series: 6-1C. Type 122540.
11. Medical diagnostic (examination) disposable non-sterile gloves “BENOVY” natural latex textured High Density powder free with inner and outer chlorination (beaded cuff), color: natural, size: XS, S, M, L, XL. Series: 85-2C. Type 122540.

- Medical diagnostic (examination) disposable non-sterile gloves "BENOVY" natural latex textured High Density powder free with inner and outer chlorination (extended beaded cuff), color: black, size: XS, S, M, L, XL. Series: 95-2C. Type 122540.
13. Medical diagnostic (examination) disposable non-sterile gloves "BENOVY" natural latex textured High Density powder free with inner and outer chlorination (beaded cuff), color: natural, size: XS, S, M, L, XL. Series: 10-2C. Type 122540.
14. Medical diagnostic (examination) disposable non-sterile gloves "BENOVY" natural latex textured Ultra Density powder free with inner and outer chlorination (extended beaded cuff), color: black, dark blue, size: XS, S, M, L, XL. Series: 20-2C. Type 122540.

1.2 Manufacturer (producer) of the medical device

Manufacturer (producer) of the medical device:

Hycare International Co., Ltd.

1197, Moo 3, Asia Highway, Khuanlang, Hat Yai Songkhla, 90110, Thailand

Address of the medical device manufacture:

1197, Moo 3, Asia Highway, Khuanlang, Hat Yai Songkhla, 90110, Thailand

1.3 Purpose of the medical device

The gloves are intended for single use as protective barrier on healthcare professional's hands during examination/treatment of patient or for other medical goals.

1.4 Scope of the medical device

Gloves are used during medical researches, diagnostic, therapeutic, dental and treatment procedures, sanitation and laboratory researches, as well as handling contaminated medical materials/devices/surfaces, providing hand hygiene, as also protection of patient and healthcare professional from cross-contamination.

1.5 Intended users of the medical device

Medical personnel, who carry out medical researches, diagnostic, therapeutic, dental and treatment procedures, sanitation and laboratory researches, as well as handling contaminated medical materials/devices/surfaces and other activities, which require compliance with the requirements for protection of patient/healthcare professional from cross-infection.

1.6 Application conditions of the medical device

Medical examination disposable gloves are used in medical or healthcare facilities indoors in a dry place, away from direct sunlight, at temperature from 0°C to 25°C and relative humidity from 65% up to 85% in accordance with the instruction for use on medical device.

Gloves for external use only (the device is non-invasive).

1.7 Type of contact with human body

Short-time surface contact with intact skin and mucous membranes of the human body.

1.8 Principle of operation of the medical device

The base of examination natural latex disposable gloves restricts the movement of infectious substances, increases the effectiveness of providing antibacterial protection in a medical facility by protecting a medical worker from possible infections when examining / treating a patient and during other sanitary actions, as well as contributes to the reliable capture of medical instruments and increases the safety of work.

Due to the differences in the method of surface treatment, its texture, glove thickness, cuff length, and additional additives in the latex, the medical device allows you to capture almost all areas of diagnostics and therapy, including dentistry, medical, sanitary, laboratory and other medical procedures and researches.

Gloves with a smooth surface provide high tactile sensitivity, and a powdering agent makes it easy to put on and take off the gloves, as well as prevents them from sticking together during storage.

Variants of non-powdered medical products with an internal polymer coating are characterized by the fact that such coating gives additional strength to the glove, isolates the skin of a medical worker from effects of natural latex proteins and significantly eases putting on and removing, as well as prevents the gloves from sticking together.

Non-powdered gloves with internal and external polymer coating provide a protective barrier thus protecting the skin of both the medical worker and the patient from effects of natural latex proteins, which contributes to the lowest risk of allergic reactions.

Due to the natural properties of Aloe Vera, versions of medical devices with the addition of this extract help to moisturize the hands skin.

Non-powdered chlorinated gloves are treated with chlorine on the inside, which makes it possible to smooth the surface of the medical device, making it easy to put on and take off the gloves, and it also prevents them from sticking together during storage.

Non-powdered gloves with internal and external chlorination are treated with chlorine on both sides, smoothing them out, making it easy to put on and take off the gloves, and preventing them from sticking together during storage.

High-density gloves are characterized by an increased wall thickness, contributing to increased resistance to damage, and are used in conditions of high infection risk during procedures, where a possibility of mechanical damage to the gloves exists, and more reliable hand protection is required.

Versions of the medical device with a textured surface ensure reliable grip and retention of medical instruments.

Gloves with an extended cuff provide improved protection of the wrist and forearm from possible contamination with infectious medical materials.

1.9 Indications for use of the medical device

The indications for use of the medical device are the necessity of conduction of medical researches, diagnostic, therapeutic, dental and treatment procedures, sanitation and laboratory researches, as well as handling contaminated medical materials/devices/surfaces.

1.10 Contradictions of the medical device

Contradictions of the medical device are as follows:

1. The patient or healthcare professional has an allergic reaction to natural latex proteins;
2. The patient or healthcare professional has an allergic reaction to the powder (applicable only to powdered gloves);
3. The patient or healthcare professional has an allergic reaction to residual process chemicals in the medical device.

1.11 Shelf life of the medical device

The established shelf life of the examination disposable non-sterile powdered natural latex gloves is 5 years.

The established shelf life of the examination disposable non-sterile powder free natural latex gloves is 3 years.

2 Warning and precautions

- The device is intended for single use;
- For individual use only;
- Every time it is necessary to perform hygienic treatment of your hands before putting on or changing gloves, as well as after removing the gloves;
- It is necessary to choose the gloves of the right size – they should be easily put on and should not impede movement or squeeze the hands. Do not use in case of size mismatch;
- It is necessary to choose the gloves depending on the type of the performed manipulations;
- It is recommended to touch only those areas of the glove that are designed for wrists (top edge of the cuff);
- It is necessary to change the gloves in case of any defect, puncture, damage, contamination, fluid contact with skin under the glove, as well as stickiness;
- If one glove is damaged, change both ones;
- Depending on the type of the performed manipulation, to keep an optimal level of protection, it is recommended to change gloves every 60-120 minutes;
- Do not wet and keep away from moisture;
- It is necessary to put the gloves on dry hands only;
- Avoid exposure to direct sunlight;
- The composition of the medical device includes natural rubber latex that may cause an allergic reaction;
- Powder may cause an allergic reaction (applicable only to powdered gloves);
- Gloves contain residual process chemicals that may cause an allergic reaction;
- Gloves do not guarantee 100% hand protection against punctures, cuts and damage caused by exposure to chemicals.
- It is recommended to consult healthcare professional in the event of dermal hypersensitivity reactions, when using the device.

3 Instruction for use of the medical device

Preliminary hygienic treatment of your hands is required before use of the gloves (treat hands with skin antiseptic or wash with soap and water). Gloves are put on dry hands only.

How to put the gloves on:

- 1) Take a glove from the box. It is recommended to touch the bounded surface of the glove corresponding to the wrist (upper edge of the cuff).
- 2) Pull on the first glove.
- 3) Take the second glove with your bare hand, touching the bounded surface of the glove corresponding to the wrist.
- 4) To avoid touching the skin of the forearm with a gloved hand, rotate the outer surface of the glove to be put on with folded fingers in the glove so that you could pull the glove over the other hand.
- 5) Having put the gloves on, it is required to touch only those objects provided for by the indications for use and conditions of the use of the gloves.
- 6) To avoid contamination, it is required to change gloves before proceeding with the new procedure.

The used gloves should be put off carefully, avoiding splashing from the glove's surfaces where microbial contamination of hands and environment is possible.

How to take the gloves off:

- 1) To take a glove off, pick it up at the level of the wrist without touching the skin of the forearm and pull it off your hand, turning it inside out.
- 2) Holding the removed glove with the gloved hand, it is necessary to place the fingers of your bare hand between the glove and wrist. Remove the second glove by rolling it down the hand and put into the first glove.
- 3) Discard the gloves and wash your hands thoroughly.

4 Transportation, storage and disposal conditions of the medical device

4.1 Transportation conditions of the medical device

The transportation of the device may be made by all means of transport (sea, land, railway or air transport) according to the rules of transportation, which are applied to such means of transport. During transportation, the storage conditions should not be violated.

4.2 Storage conditions of the medical device

Devices should be stored in a dry place protected from direct sunlight, heat, ignition and ultraviolet sources, as well as ozone, at temperature from 0°C to 25°C and relative humidity – from 65% up to 85%.

Products should be handled with caution when they are stored at temperatures below 15°C, as during careless handling, they can be stiff and sensitive to deformation. After storage at a low temperature, gloves should be warmed up to 30°C before use throughout all product volume.

4.3 Disposal conditions of the medical device

Disposal of the medical device in the Russian Federation is carried out by the medical organization as class B waste in accordance with SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of urban and rural settlements, for water facilities, drinking water and drinking water supply, atmospheric air, soils, living quarters, the operation of industrial, public premises, the organization and implementation of sanitary and anti-epidemic (preventive) measures" and approved guideline of the organization.

5 Cleaning and disinfection methods of the medical device

Examination gloves are disposable medical device. Reuse of the medical device is a gross violation of application methods.

6 Specification of the medical device

6.1 Description of the medical device

Form of the medical device. Medical examination gloves have straight fingers and are ambidextrous (non-anatomic). The glove's cuff may be extended or has a standard length. Rolled cuff.

Surface of the medical device.

1. Powdered gloves - modified corn starch is used.
2. Powder-free gloves:
 - 2.1. Chlorinated (treated with chlorine on the inside wall of gloves);
 - 2.2. With internal and external chlorination;
 - 2.3. With internal polymer coating;
 - 2.4. With internal and external polymer coating.

The texture of the external surface of the medical device. Gloves may have a smooth or textured surface.

Sizes of the medical device. Depending on the size, table 1 gives the values of width, length and thickness of the medical device.

Table 1 – Sizes of the medical device:

Value name	Size				
	XS	S	M	L	XL
Width W (figure 1), mm	≤ 80	80±10	95±10	110±10	≥110
Length l (figure 1), mm, min	220		230		
Thickness (in points, shown on figure 2), mm, min: - smooth surface - textured surface	0.08 0.11				
Thickness of rolled cuff, mm, max	2.50				


Figure 1 – Points for width and length measurement of the glove


Figure 2 – Points for wall thickness measurement of the glove

6.2 Physical specification of the medical device

6.2.1 Physical and mechanical performance of the medical device

Physical and mechanical parameters of the medical device are represented in table 2.

Table 2 – Values of physical and mechanical parameters of gloves:

DESCRIPTION	SPECIFICATION
<u>1. Before aging:</u>	
- conditional strength, MPa, min	18
- elongation at break, %, min	650
- force at break, N, min	7.0
<u>2. After aging:</u>	
- conditional strength, MPa, min	14
- elongation at break, %, min	500
- force at break, N, min	6.0

6.2.2 Inspection levels and acceptable quality level of the medical device

For the purposes of quality control medical gloves are selected and examined in accordance with standard ISO 2859-1.

Hermeticity. Medical gloves are hermetic. Acceptable quality level (AQL) is ≤ 1.5 .

Powder. The content of residual powder is determined in accordance with the requirements of ASTM D 6124.

- For powdered gloves, the value obtained does not exceed 10 mg/dm² per glove.
- For powder-free gloves, the value obtained does not exceed 2 mg per glove.

The protein content in latex gloves. The content of water extractable proteins is not more than 200 µg/dm² when tested in accordance with ASTM D 5712.

7 Adverse reactions

Manifestation of dermal hypersensitivity reactions, when using the device is possible.

8 Circumstances, in which the consumer must consult healthcare professional

The patient is recommended to consult healthcare professional in the event of dermal hypersensitivity reactions, when using the device.

9 Manufacturer's warranties

The manufacturer (producer) guarantees preservation of all the parameters and characteristics of the medical device taking into account acceptable quality level (AQL), as well as the quality of the medical gloves within all shelf life from the date of manufacture upon the condition of package integrity, observance of transportation, storage and operation conditions in accordance with the instruction.

10 Data on Authorized representative of the manufacturer

Please contact the Authorized representative of the manufacturer for any and all purposes, related to circulation of the medical device on the territory of the Russian Federation:

“BENOVY” Company Limited (“BENOVY” Co., Ltd.)

Russia, 194292, St. Petersburg, 1, 5-y Verkhniy pereulok, litera J, office 205.

tel. +7 (812) 309-16-32

www.benovy-company.ru

e-mail: benovy@benovy-company.ru

11 Scope of delivery of the medical device

Gloves can be delivered in the following scopes of supply:

- bulk packaging – 2000 gloves, 1000 pairs in consumer package;
- 400 gloves, i.e. 200 pairs in consumer package;
- 200 gloves, i.e. 100 pairs in consumer package;
- 100 gloves, i.e. 50 pairs in consumer package;
- 50 gloves, i.e. 25 pairs in consumer package;
- 40 gloves, i.e. 20 pairs in consumer package;
- 30 gloves, i.e. 15 pairs in consumer package;
- 20 gloves, i.e. 10 pairs in consumer package;
- 10 gloves, i.e. 5 pairs in consumer package;
- 4 gloves, i.e. 2 pairs in consumer package.

12 List of national standards, applied by manufacturer (producer) of the medical device

Table 3 - List of standards, applied by manufacturer:

Designation	Name
ISO 9001:2015	Quality management systems. Requirements.
ISO 13485:2016	Medical devices. Quality management systems. Requirements for regulatory purposes.
ISO 14971	Medical devices – Application of risk management to medical devices.
ISO 2859-1	Statistical methods. Sampling procedures for inspection by attributes. Part 1. Sampling schemes indexed by acceptance quality limit for lot-by-lot inspection.
ISO 15223-1:2012	Medical devices. Symbols to be used with medical device labels, labelling, and information to be supplied. Part 1. General requirements.
ISO 11193-1	Single-use medical examination gloves. Part 1: Specification for gloves made from rubber latex or rubber solution.
ASTM D 3578	Standard specification for rubber examination gloves.
ASTM D 6124	Standard test method for residual powder on medical gloves.
ASTM D 5712	Standard test method for analysis of aqueous extractable protein in latex, natural rubber, and elastomeric products using the modified Lowry method.

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