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Note: Model 28601 is the blade most commonly used for left internal mammary artery (LIMA) harvest. Model 28605 is the blade most commonly used for general thoracotomy retraction. All other blades were designed to accommodate various patient anatomies.

Note: Models 28604B and 28606B are inferior blades designed for women and obese patients. Use these inferior blades with the superior blades (Models 28604 and 28606).

2 Indications for use

The ThoraTrak minimally invasive cardiac surgery (MICS) retractor system is intended to provide surgical access for minimally invasive cardiothoracic procedures, including minimally invasive coronary artery bypass grafting (CABG) surgery and LIMA harvest, by retraction of soft and bony tissue.

3 Precautions

Proper surgical procedures and techniques are the responsibility of the medical professional. Surgeons must evaluate the appropriateness of any procedure based on their own medical training and experience.

4 Adverse effects

None known.

5 Processing and reprocessing

Medtronic has validated the following instructions for preparing this medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials, and personnel in the processing facility, achieves the desired result. Processing requires validation and routine monitoring. Likewise, any deviation by the processor from the following instructions should be properly evaluated for effectiveness and potential adverse consequences. Cleaning may be done manually, or it may be automated according to these Instructions for Use or an equivalent, validated method.

5.1 Cleaning

Warning: Thoroughly clean the device to eliminate soil prior to sterilization.

Warning: Do not use solvent-based cleaners (for example, acetone or toluene).

- Do not allow contaminated devices to dry before implementing any cleaning procedure.
- Reprocess devices within 2 hours after use. If transport to the processing area will be delayed, place the device in a covered container with an enzymatic detergent to prevent drying.
- If the device contains multiple components, completely disassemble the components prior to cleaning per these Instructions for Use.
- Examine the device for cracking or degradation prior to use. Do not use the device if any components show signs of cracking or degradation.
- Thermal disinfection is not required because the devices are terminally sterilized.

5.1.1 Manual cleaning

Medtronic has established the following manual cleaning procedure in Table 1. These steps must be performed within a maximum of 2 hours after use.

Note: If the device has movable parts, ensure that all surface areas have been thoroughly cleaned.

Table 1. Manual cleaning instructions

Step	Process	Temperature	Cleaning Instructions
1	Rinsing	>27°C (>81°F)	Remove gross soil using running water for 1 min. Use a soft-bristle brush (for example, a nylon toothbrush) to clean the device thoroughly.
2	Soaking	>28°C (>82°F)	Submerge the device completely for a minimum of 5 min with 15.6 mL/1 L (or 2 ounces/1 gallon) of an enzymatic detergent (for example, ENZOL™) and water. Remove visible soil with a soft-bristled brush.
3	Ultrasonic clean	>27°C (>81°F)	Submerge the device completely in an ultrasonic cleaner (for example, Branson™) with 7.8 mL/1 L (or 1 ounce/1 gallon) of an enzymatic detergent (for example, ENZOL) and water. Sonicate for 10 min.
4	Rinsing	>27°C (>81°F)	Rinse the device with running water for 1 min.
5	Drying	N/A	Dry with a clean, lint-free wipe.
6	Inspection	N/A	Visually inspect each device for any remaining soil or moisture. If soil remains, repeat the process.

5.1.2 Automated cleaning

Medtronic has established the following automated cleaning procedure. The steps in Table 2 must be performed within a maximum of 2 hours after use.

Note: If the device or component has movable parts, ensure that all surface areas have been thoroughly cleaned.

Table 2. Pretreatment instructions

Step	Process	Water temperature	Cleaning instructions
1	Rinsing	>25°C (>77°F)	Remove gross soil using running tap water.
2	Soaking	>25°C (>77°F)	Submerge the device completely for a minimum of 1 min with 7.8 mL/1 L (or 1 ounce/1 gallon) of an enzymatic detergent (for example, ENZOL) and water. Remove visible soil with a soft-bristled brush (for example, a nylon toothbrush).
3	Rinsing	>25°C (>77°F)	Rinse the device for a minimum of 2 min.
4	Inspection	N/A	Visually inspect the device for any remaining soil. If needed, repeat the process.

Medtronic used the Polystica™ family of cleaners according to the manufacturer's recommendations to validate the automated cleaning process. It is the responsibility of the processor to ensure that the processing is conducted in accordance with a validated method. Any deviation by the processor from these recommendations should be evaluated.

1. Place the device in an automated washer (for example, Steris Reliance™ Genfore™ Washer/Disinfector).

Note: Avoid contact between devices while loading the washer.

2. Set the automated washer to run the parameters identified in Table 3 and allow the automated washer to complete the wash cycle.

Table 3. Automated washer cleaning parameters

Treatment	Time (min)	Temperature	Cleaning agent
Enzymatic wash	04:00	≥60°C (≥140°F)	Steris Polystica™ Ultra Concentrate Enzymatic Cleaner diluted as 1.0 mL/1 L (or 0.125 ounce/1 gallon).
Wash	02:00	≥60°C (≥140°F)	Steris Polystica™ Ultra Concentrate Neutral Detergent diluted as 1.0 mL/1 L (or 0.125 ounce/1 gallon).
Rinse	02:00	≥60°C (≥140°F)	N/A
Dry	15:00	≥82°C (≥180°F)	N/A
Inspection	N/A	N/A	Visually inspect each device for any remaining soil or moisture. If soil remains, repeat the process. If needed, dry the device with filtered, compressed air or a lint-free wipe.

5.2 Steam sterilization and resterilization

- Use standard packaging material. Ensure that the pack is large enough to contain the device without stressing the seals. When validating the sterilization processes, Medtronic used the appropriate accessory tray for each device.

- Medtronic validated steam cycles when the devices were wrapped by CSR (Central Supply Room) wrap. However, the steam sterilization process can be completed when the devices are either unwrapped or wrapped in steam sterilization compatible materials.
- When sterilizing multiple devices in 1 autoclave cycle, do not exceed the sterilizer's maximum load.
- Examine the device for cracking or degradation prior to use. Do not use the device if components show signs of cracking or degradation.
- Some non-US health care authorities recommend sterilization cycle parameters that minimize the potential risk of transmitting Creutzfeldt-Jakob Disease (CJD). This recommendation is especially important for surgical instruments that could come into contact with the central nervous system.

Table 4. Sterilization cycle parameters

Cycle type	Prevacuum (dynamic-air-removal)	Prevacuum (dynamic-air-removal)	Prevacuum (dynamic-air-removal) for CJD ^a
Temperature	132°C (270°F)	135°C (275°F)	134°C (273°F)
Exposure time	4 min	3 min	18 min
Dry time ^b	20 min	16 min	30 min

^a Medtronic recommends incinerating devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE)/CJD diagnosis. *NHS Estates HTM 2010 Parts 4 & 6: Appendix 2, Items contaminated with TSE Agents* and *WHO Infection Control Guidelines for Transmissible Spongiform Encephalopathies* refer to a TSE decontamination cycle using a steam autoclave at a temperature of 134°C to 137°C (273°F to 279°F) for a single cycle of 18 min or multiple cycles totaling 18 min (for example, six 3 min cycles).

^b The minimum dry times were validated using sterilizers with vacuum drying capabilities. Drying cycles using ambient atmospheric pressure may require longer dry times. Refer to the sterilizer manufacturer's recommendations.

6 Instructions for use

Note: The blades are marked with a number that corresponds to the last digit of the model number. For example, 28603 is marked "3"; 28606B is marked "6B." Each blade is marked with a unique serial number.

6.1 Blade attachment

1. Attach the blades to the retractor rack by sliding the blade mount into the retractor rack arm.
2. Remove the blades by sliding the blade mount out of the retractor rack arm.

6.2 LIMA harvest

1. Place the retractor, with LIMA blades attached, into the thoracotomy. The large LIMA blade with the hole at the end should be oriented cephalad to the smaller blade.
2. Turn the retractor knob to open the retractor and spread the thoracotomy.
3. Apply upward tension on the large LIMA blade at the hole at the end of the blade to tilt the retractor upward. This action provides access to the LIMA and enables harvesting.

Note: The retractor is designed to work with the Rultract™ retractor system to apply upward tension on the large LIMA blade. Other means or methods for applying upward tension may be used.

4. Once the LIMA harvest is complete, release upward tension from the large LIMA blade.
5. Depress the release lever to close the retractor.
6. Remove the retractor from the thoracotomy.

6.3 General thoracotomy retraction

1. Attach thoracotomy blades to the retractor.
2. Place the retractor into the thoracotomy.
3. Turn the retractor knob to open the retractor and spread the thoracotomy.
4. Upon completion of the procedure, depress the release lever to close the retractor.
5. Remove the retractor from the thoracotomy.

7 Disclaimer of warranty (global)

ALTHOUGH THE PRODUCT HAS BEEN DESIGNED, MANUFACTURED, AND TESTED PRIOR TO SALE UNDER CAREFULLY CONTROLLED CONDITIONS, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION FOR A VARIETY OF REASONS. MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE, OR MALFUNCTION OF THE PRODUCT, WHETHER THE CLAIM IS BASED UPON WARRANTY, CONTRACT, TORT, OR OTHERWISE. EXCEPT AS OTHERWISE AGREED TO IN WRITING BY AUTHORIZED MEDTRONIC PERSONNEL, NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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© 2016 Medtronic
The following list includes trademarks or registered trademarks of Medtronic in the United States and possibly in other countries. All other trademarks are the property of their respective owners.
(Kluge Design, Inc.)