APPROVE

Poly Medicure Limited

Sr. GM-QUALITY

(position)

SS RAWAT

(name)

SS Rawat

(signature)

Stamp

INSTRUCTIONS FOR USE

**Spinal Quincke and Pencil Point needle for anesthesia in various designs
manufactured by
Poly Medicure Limited**

24 FEB 2020

ATTESTED

W
**NOTARY PUBLIC
(INDIA)**

2020

1. NAME AND COMPOSITION OF THE MEDICAL DEVICE

Spinal Quincke and Pencil Point needles for anesthesia in various designs:

I. Spinal Quincke needle for anesthesia in various designs:

1. Spinal Quincke needle for anesthesia – size 16G (1.6 mm), length 89 mm in the composition:

1.1. Spinal Quincke needle for anesthesia - size 16G (1.6 mm), length 89 mm

1.2. Instructions for use.

2. Spinal Quincke needle for anesthesia - size 17G (1.4 mm), length 89 mm in the composition:

2.1 Spinal Quincke needle for anesthesia - size 17G (1.4 mm), length 89 mm

2.2 Instructions for use.

3. Spinal Quincke needle for anesthesia - size 18G (1.2 mm), length 89 mm in the composition:

3.1 Spinal Quincke needle for anesthesia - size 18G (1.2 mm), length 89 mm

3.2 Instructions for use.

4. Spinal Quincke needle for anesthesia - size 19G (1.1 mm), length 89 mm in the composition:

4.1 Spinal Quincke needle for anesthesia - size 19G (1.1 mm), length 89 mm

4.2 Instructions for use.

5. Spinal Quincke needle for anesthesia - size 20G (0.9 mm), length 89 mm in the composition:

5.1 Spinal Quincke needle for anesthesia - size 20G (0.9 mm), length 89 mm

5.2 Instructions for use.

6. Spinal Quincke needle for anesthesia - size 21G (0.8 mm), length 89 mm in the composition:

6.1 Spinal Quincke needle for anesthesia - size 21G (0.8 mm), length 89 mm

6.2 Instructions for use.

7. Spinal Quincke needle for anesthesia - size 22G (0.7 mm), length 89 mm in the composition:

7.1 Spinal Quincke needle for anesthesia - size 22G (0.7 mm), length 89 mm

7.2 Instructions for use.

8. Spinal Quincke needle for anesthesia - size 23G (0.6 mm), length 89 mm in the composition:

8.1 Spinal Quincke needle for anesthesia - size 23G (0.6 mm), length 89 mm

8.2 Instructions for use.

9. Spinal Quincke needle for anesthesia - size 24G (0.55 mm), length 89 mm in the composition:

9.1 Spinal Quincke needle for anesthesia - size 24G (0.55 mm), length 89 mm

9.2 Instructions for use.

10. Spinal Quincke needle for anesthesia - size 25G (0.5 mm), length 89 mm in the composition:

10.1 Spinal Quincke needle for anesthesia - size 25G (0.5 mm), length 89 mm

10.2 Instructions for use.

11. Spinal Quincke needle for anesthesia - size 26G (0.45 mm), length 89 mm in the composition:


11.1 Spinal Quincke needle for anesthesia - size 26G (0.45 mm), length 89 mm

11.2 Instructions for use.

12. Spinal Quincke needle for anesthesia - size 27G (0.4 mm), length 89 mm in the composition:

12.1 Spinal Quincke needle for anesthesia - size 27G (0.4 mm), length 89 mm

12.2 Instructions for use.

ATTESTED

(Signature)

NOTARY PUBLIC
(INDIA)

24 FEB 2020

33675
21.05.2020

13. Spinal Quincke needle for anesthesia - size 28G (0.36 mm), length 89 mm in the composition:

- 13.1 Spinal Quincke needle for anesthesia - size 28G (0.36 mm), length 89 mm
- 13.2 Instructions for use.

14. Spinal Quincke needle for anesthesia - size 29G (0.33 mm), length 89 mm in the composition:

- 14.1 Spinal Quincke needle for anesthesia - size 29G (0.33 mm), length 89 mm
- 14.2 Instructions for use.

15. Spinal Quincke needle for anesthesia - size 30G (0.3 mm), length 89 mm in the composition:

- 15.1 Spinal Quincke needle for anesthesia - size 30G (0.3 mm), length 89 mm
- 15.2 Instructions for use.

16. Spinal Quincke needle for anesthesia - size 22G (0.7 mm), length 89 mm with introducer needle, size 18G (1.2 mm), length 38 mm in the composition:

- 16.1 Spinal Quincke needle for anesthesia - size 22G (0.7 mm), length 89 mm
- 16.2 Introducer needle, size 18G (1.2 mm), length 38.
- 16.3 Instructions for use.

17. Spinal Quincke needle for anesthesia - size 23G (0.6 mm), length 89 mm with introducer needle, size 19G (1.1 mm), length 38 mm in the composition:

- 17.1 Spinal Quincke needle for anesthesia - size 23G (0.6 mm), length 89 mm
- 17.2 Introducer needle, size 19G (1.1 mm), length 38.
- 17.3 Instructions for use.

18. Spinal Quincke needle for anesthesia - size 24G (0.55 mm), length 89 mm with introducer needle, size 20G (0.9 mm), length 38 mm in the composition:

- 18.1 Spinal Quincke needle for anesthesia - size 24G (0.55 mm), length 89 mm
- 18.2 Introducer needle, size 20G (0.9 mm), length 38.
- 18.3 Instructions for use.

19. Spinal Quincke needle for anesthesia - size 25G (0.5 mm), length 89 mm with introducer needle, size 21G (0.8 mm), length 38 mm in the composition:

- 19.1 Spinal Quincke needle for anesthesia - size 25G (0.5 mm), length 89 mm
- 19.2 Introducer needle, size 21G (0.8 mm), length 38.
- 19.3 Instructions for use.

20. Spinal Quincke needle for anesthesia - size 26G (0.45 mm), length 89 mm with introducer needle, size 21G (0.8 mm), length 38 mm in the composition:

- 20.1 Spinal Quincke needle for anesthesia - size 26G (0.45 mm), length 89 mm
- 20.2 Introducer needle, size 21G (0.8 mm), length 38.
- 20.3 Instructions for use.

21. Spinal Quincke needle for anesthesia - size 27G (0.4 mm), length 89 mm with introducer needle, size 22G (0.7 mm), length 38 mm in the composition:

- 21.1 Spinal Quincke needle for anesthesia - size 27G (0.4 mm), length 89 mm
- 21.2 Introducer needle, size 22G (0.7 mm), length 38.
- 21.3 Instructions for use.

Spinal Pencil Point needle for anesthesia in various designs:

1. Spinal Pencil Point needle for anesthesia - size 16G (1.6 mm), length 89 mm in the composition:

- 1.1. Spinal Pencil Point needle for anesthesia - size 16G (1.6 mm), length 89 mm.
- 1.2 Instructions for use.

2. Spinal Pencil Point needle for anesthesia - size 17G (1.4 mm), length 89 mm in the composition:

- 2.1. Spinal Pencil Point needle for anesthesia - size 17G (1.4 mm), length 89 mm.
- 2.2 Instructions for use.

3. Spinal Pencil Point needle for anesthesia - size 18G (1.2 mm), length 89 mm in the

24 FEB 2020

ATTESTED

NOTARY PUBLIC
(INDIA)

composition:

3.1. Spinal Pencil Point needle for anesthesia – size 18G (1.2 mm), length 89 mm.

3.2 Instructions for use.

4. Spinal Pencil Point needle for anesthesia – size 19G (1.1 mm), length 89 mm in the composition:

4.1. Spinal Pencil Point needle for anesthesia – size 19G (1.1 mm), length 89 mm.

4.2 Instructions for use.

5. Spinal Pencil Point needle for anesthesia – size 20G (0.9 mm), length 89 mm in the composition:

5.1. Spinal Pencil Point needle for anesthesia – size 20G (0.9 mm), length 89 mm.

5.2 Instructions for use.

6. Spinal Pencil Point needle for anesthesia – size 21G (0.8 mm), length 89 mm in the composition:

6.1. Spinal Pencil Point needle for anesthesia – size 21G (0.8 mm), length 89 mm.

6.2 Instructions for use.

7. Spinal Pencil Point needle for anesthesia – size 22G (0.7 mm), length 89 mm in the composition:

7.1. Spinal Pencil Point needle for anesthesia – size 22G (0.7 mm), length 89 mm.

7.2 Instructions for use.

8. Spinal Pencil Point needle for anesthesia – size 23G (0.6 mm), length 89 mm in the composition:

8.1. Spinal Pencil Point needle for anesthesia – size 23G (0.6 mm), length 89 mm.

8.2 Instructions for use.

9. Spinal Pencil Point needle for anesthesia – size 24G (0.55 mm), length 89 mm in the composition:

9.1. Spinal Pencil Point needle for anesthesia – size 24G (0.55 mm), length 89 mm.

9.2 Instructions for use.

10. Spinal Pencil Point needle for anesthesia – size 25G (0.5 mm), length 89 mm in the composition:

10.1. Spinal Pencil Point needle for anesthesia – size 25G (0.5 mm), length 89 mm.

10.2 Instructions for use.

11. Spinal Pencil Point needle for anesthesia – size 26G (0.45 mm), length 89 mm in the composition:

11.1. Spinal Pencil Point needle for anesthesia – size 26G (0.45 mm), length 89 mm.

11.2 Instructions for use.

12. Spinal Pencil Point needle for anesthesia – size 27G (0.4 mm), length 89 mm in the composition:

12.1. Spinal Pencil Point needle for anesthesia – size 27G (0.4 mm), length 89 mm.

12.2 Instructions for use.

13. Spinal Pencil Point needle for anesthesia – size 28G (0.36 mm), length 89 mm in the composition:

13.1. Spinal Pencil Point needle for anesthesia – size 28G (0.36 mm), length 89 mm.

13.2 Instructions for use.

14. Spinal Pencil Point needle for anesthesia – size 29G (0.33 mm), length 89 mm in the composition:

14.1. Spinal Pencil Point needle for anesthesia – size 29G (0.33 mm), length 89 mm.



14.2 Instructions for use.

15. Spinal Pencil Point needle for anesthesia – size 30G (0.3 mm), length 89 mm in the composition:

15.1. Spinal Pencil Point needle for anesthesia – size 30G (0.3 mm), length 89 mm.

15.2 Instructions for use.

16. Spinal Pencil Point needle for anesthesia – 22G (0.7 mm), length 89 mm with intro-

TESTED

(w)

NOTARY PUBLIC (INDIA)

24 FEB 2020

ducer needle, size 18G (1.2 mm), length 38 mm in the composition:

16.1. Spinal Pencil Point needle for anesthesia – 22G (0.7 mm), length 89 mm.

16.2. Introducer needle, size 18G (1.2 mm), length 38 mm.

16.3 Instructions for use.

17. Spinal Pencil Point needle for anesthesia – 23G (0.6 mm), length 89 mm with introducer needle, size 19G (1.1 mm), length 38 mm in the composition:

17.1. Spinal Pencil Point needle for anesthesia – 23G (0.6 mm), length 89 mm.

17.2. Introducer needle, size 19G (1.1 mm), length 38 mm.

17.3 Instructions for use.

18. Spinal Pencil Point needle for anesthesia – 24G (0.55 mm), length 89 mm with introducer needle, size 20G (0.9 mm), length 38 mm in the composition:

18.1. Spinal Pencil Point needle for anesthesia – 24G (0.55 mm), length 89 mm.

18.2. Introducer needle, size 20G (0.9 mm), length 38 mm.

18.3 Instructions for use.

19. Spinal Pencil Point needle for anesthesia – 25G (0.5 mm), length 89 mm with introducer needle, size 21G (0.8 mm), length 38 mm in the composition:

19.1. Spinal Pencil Point needle for anesthesia – 25G (0.5 mm), length 89 mm.

19.2. Introducer needle, size 21G (0.8 mm), length 38 mm.

19.3 Instructions for use.

20. Spinal Pencil Point needle for anesthesia – 26G (0.45 mm), length 89 mm with introducer needle, size 21G (0.8 mm), length 38 mm in the composition:

20.1. Spinal Pencil Point needle for anesthesia – 26G (0.45 mm), length 89 mm.

20.2. Introducer needle, size 21G (0.8 mm), length 38 mm.

20.3 Instructions for use.

21. Spinal Pencil Point needle for anesthesia – 27G (0.4 mm), length 89 mm with introducer needle, size 22G (0.7 mm), length 38 mm in the composition:

21.1. Spinal Pencil Point needle for anesthesia – 27G (0.4 mm), length 89 mm.

21.2. Introducer needle, size 22G (0.7 mm), length 38 mm.

21.3 Instructions for use.

2. MANUFACTURER OF THE MEDICAL DEVICE

Poly Medicure Limited (Medcomp), 232-B, 3rd Floor, Okhla Industrial Estate, Phase III, New Delhi – 110 020 India

Tel/Fax: + 91- 11- 33550700/+ 91-11-26321894

info@polymedicure.com

Production site of medical device:

Plot № 115-116, Sector 59, HSIIDC, Industrial Area, Ballabhgarh, Faridabad-121004, Haryana, India.

3. INTENDED USE OF THE MEDICAL DEVICE DETERMINED BY THE MANUFACTURER

Spinal needles in various designs (hereinafter – infusion lines) are used for lumbar puncture through the dura mater of the spinal column to administer anesthetics used for general anesthesia of the human body

Indications

- surgeries below the navel level;
- gynecological and urological surgeries;

ATTESTED

NOTARY PUBLIC
(INDIA)

24 FEB 2020

- surgical manipulations on the lower extremities, for example, the treatment of varicose veins;
- perineal surgery;
- pain relief during childbirth and cesarean section;
- alternative to general anesthesia in case of contraindications (old age, somatic pathology, allergy to anesthetics).

Absolute Contraindications

- patient's refusal;
- blood diseases accompanied by increased bleeding, the use of anticoagulants prior to the surgery (high risk of bleeding);
- inflammatory lesions of the skin at the expected puncture site;
- severe condition of the patient (shock, acute blood loss, cardiovascular, pulmonary failure, sepsis);
- allergy to local anesthetics used for analgesia;
- infectious diseases of the nervous system (meningitis, arachnoiditis, encephalitis, myelitis);
- intracranial hypertension;
- exacerbation of herpes virus infection;
- severe degree of cardiac arrhythmias and blockades.

Relative Contraindications

- spinal column deformity, which increases the risk of complications and makes anesthesia dangerous to life and health;
- predicted extensive blood loss during future surgery;
- severe degree of fetal distress when choosing a method of delivery;
- signs of infectious disease, fever;
- some diseases of the NS (epilepsy, radiculitis with radicular syndrome, vascular lesions of the brain, polio, chronic headache, multiple sclerosis);
- patient's emotional instability, mental disorders (persons unable to lie still during surgery);
- stenosis of the aortic valve of the heart;
- therapy with acetylsalicylic acid and other disaggregants (risk of bleeding);
- history of spinal cord injuries;
- possible expansion of the scope and time of the surgery, for example, surgical removal of tumors, when the surgeon's tactics may change depending on the picture on the operating table;
- child age.

Risks of medical device, contraindications, expected and predictable side effects associated with the intended use of medical device:

1. Before using a medical device, please, carefully read the instructions for use.
2. Do not use the product if it is damaged, of inadequate quality or the packaging is opened.
3. Violation of the storage conditions may impair performance of the product.
4. Disposable product; repeated sterilization may cause malfunction.
5. Excessive relative humidity inside the sterilization chamber may affect product functionality.
6. The product is to be used only by qualified personnel, as improper use of the product may affect the patient's health.

Side effects associated with the use of the product:

ATTESTED

NOTARY PUBLIC
(INDIA)

33675
21.05.2020

1. Spinal shock can lead to the cessation of nervous activity, including loss of motor, sensory, reflex and autonomic activity.
2. Lumbar puncture with increased intracranial pressure can cause clinical disorders, for example, violation of the blood supply, causing strokes and various neurological disorders.
3. Lumbar puncture can cause bleeding, compressing the spinal cord during blood clotting, and along with skin infections can lead to complications such as cellulitis, furunculosis.
4. Postdural or lumbar puncture headache can cause severe pain, including in the back, neck and shoulders, leading to numbness of the neck.
5. Hematoma of the spinal canal can lead to accumulation of blood in the epidural region and compression of the spinal cord, which causes pain, muscle weakness, diaphragm and pelvic disorders.
6. An epidural abscess can cause fever and back pain, neurological disorders of the lower extremities and loss of diaphragm / pelvis function.
7. Allergic reactions and violations of the intended use may occur in the patients with hypersensitivity to the materials used in the product.

4. MEDICAL DEVICE CLASSIFICATION

Depending on the degree of potential risk, "Spinal needles in various designs" belongs to class **III** medical device, in accordance with European Directive 93/42/EEC.


5. DELIVERY

Spinal needles in various designs are delivered in individual package and instructions for use.

Their designs are listed in Section 1 "NAME OF THE MEDICAL DEVICE"

ATTESTED


NOTARY PUBLIC
(INDIA)

24 FEB 2020


6. INSTRUCTIONS FOR USE

Spinal anesthesia: A local anesthetic is injected into the cerebrospinal fluid, usually in the area of the lumbar spine, where it interacts with the roots of the spinal nerves and partially spinal column. The resulting anesthesia usually affects the area from the legs up to the abdomen or chest. Lumbar puncture is used to inject medications into the cerebrospinal fluid (intrathecal), in particular, to perform spinal anesthesia or chemotherapy.

The Quincke or Pencil Point needles are used for spinal anesthesia.

The incidence of post-puncture headache is higher than when using Pencil Point needles.

The frequency of post-puncture pain, since it is less traumatic to the dura matter.

Quincke needle. This needle is most often used for spinal anesthesia. The needle has a medium-length bevel with sharp edges, a sharp tip and an opening at the tip. The needle is characterized by good sharpening (cutting), the puncture with such needles is technically simpler, however, due to trauma to the dura mater, the incidence of post-puncture headache is higher than when using Pencil Point needles.

Pencil Point needle. It has a rounded, non-cutting tip and a lateral opening, in the proximal section of the needle bevel. The conical shape of the spinal needle tip without cutting edges makes it possible to puncture the solid dura mater non-traumatically. The lateral opening at the tip allows for targeted administration of an anesthetic solution. The main aim of this needle is to reduce

The instructions for use of the spinal needle with introducer needle.

To perform spinal anesthesia, the spinal needle is inserted between the spinous processes, taking into account the angle they form with respect to the axis of the spine. The spinal needle passes the following anatomical structures: skin, subcutaneous tissue, supraspinatus ligament, interspinous ligament, yellow ligament, epidural space, dura mater. These ligaments are usually very dense and even calcified in elderly patients.


In this case, it is advisable to use an additional introducer needle. The introducer is a shorter (35 - 40mm) needle of a larger diameter than the inserted spinal needle. The introducer is injected at the proposed puncture site, passes through the densest anatomical layers without deforming. Then, a spinal needle is inserted into its pavilion and, bypassing the dense tissues, a puncture of the spinal space or spinal anesthesia is performed.

Thus, the introducer is designed to pass through dense tissues and to prevent contact of the spinal needle with the skin.

Use:

- 6.1. It is suitable for spinal anesthesia.
- 6.2. Tapered needle tip provides superior tactile feedback.
- 6.3. High flow rate accelerates cerebrospinal fluid return.
- 6.4. Intensive distribution of the anesthetic immediately after injection.
- 6.5. Transparent polycarbonate needle connector provides unobstructed visualization of cerebrospinal fluid or blood.
- 6.6. Products are sterilized by gas (ethylene oxide).
- 6.7. Neither the product nor its package contains components of animal origin.

For spinal anesthesia they use 3 positions of the patient (the position of the patient when the needle is introduced):




ATTESTED

NOTARY PUBLIC
(INDIA)

24 FEB 2020

33675
21.05.2020


Needle position

24 FEB 2020

ATTESTED

NOTARY PUBLIC
(INDIA)

Lumbar puncture. (For lumbosacral access, one identifies the posterior-superior iliac spine and moves away 1 cm lower and 1 cm in medial direction. After skin anesthesia, the needle is guided through this point at an angle of 45° in the cranial and 45° medial directions and it is injected into the space between the fifth lumbar vertebra and the sacrum.)


Attention! The data given are advisory, manipulations are allowed only by trained qualified personnel.

7. TECHNICAL SPECIFICATIONS

Figure 1 and Table 1 present the sizes and tip configuration of Quincke needles, Figure 2 and Table 2 present the sizes and tip configuration of Pencil Point needles.

24 FEB 2020

ATTESTED

 NOTARY PUBLIC
 (INDIA)


Figure 1- Configuration of the Quincke needle tip

Table 1

Needle size	$A \pm 0.1$	B	C	$d_0 \pm 0.01$	$d_1 \pm 0.01$	$\alpha \pm 1^\circ$	$\theta_1, \theta_2 \pm 4^\circ$
	Edge length	Primary bevel length	Secondary bevel length	Nominal outer diameter of the needle tube	Inner diameter of the needle tube	Primary bevel angle	
16G	5.6	80 %	20%	1.60	1.20	16	35
17G	5.0			1.40	1.07		
18G	4.3			1.20	0.86		
19G	3.6			1.10	0.68		
20G	3.1			0.90	0.61		
21G	2.8			0.80	0.52		
22G	2.4			0.70	0.42		
23G	2.2			0.60	0.34		
24G	2.0			0.55	0.33		
25G	1.7			0.50	0.34		
26G	1.5			0.45	0.30		
27G	1.4			0.40	0.25		
28G	1.3			0.36	0.19		
29G	1.1			0.33	0.19		
30G	1.0			0.30	0.16		

24 FEB 2020

ATTESTED

NOTARY PUBLIC (INDIA)


Figure 2 - Configuration of the Pencil Point needle tip

Table 2

Needle size	Sharpening angle, °	Nominal outer diameter of the needle / inner diameter of the needle	Distance from the needle tip to the lateral opening	Lateral opening length
16G	30	1.60/1.20	1.20	0.6
17G	30	1.40/1.00	1.20	0.6
18G	30	1.20/0.80	1.20	0.6
19G	30	1.10/ 0.70	1.20	0.6
20G	30	0.90/0.60	1.20	0.6
21G	30	0.80/0.50	1.20	0.6
22G	30	0.70/0.40	1.20	0.6

ATTESTED

24 FEB 2020

NOTARY PUBLIC (INDIA)

33675
21.05.2020

23G	30	0.60/0.36	1.20	0.6
24G	30	0.55/0.32	1.20	0.6
25G	30	0.50/0.30	1.20	0.6
26G	30	0.45/0.28	1.20	0.6
27G	30	0.40/0.24	1.20	0.6
28G	30	0.36/0.20	1.20	0.6
29G	30	0.33/0.17	1.20	0.6
30G	30	0.30/0.14	1.20	0.6


Figure 3 – Introducer needle draft (e.g. 21G needle)

Table 3- Main technical specifications

Name	Meaning
Spinal Quincke needle - Size 16G, length 89 mm	
Size	16G (Ø1.60 mm)
Length	89 mm (±1.5mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	69N
Spinal Quincke needle - Size 17G, length 89 mm	
Size	17G (Ø1.40 mm)
Length	89 mm (±1.5mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm

24 FEB 2020

ATTESTED
NOTARY PUBLIC
(INDIA)

Needle head-tube connection stability	69N
Spinal Quincke needle - Size 18G, length 89 mm	
Size	18G (Ø1.20 mm)
Length	89 mm (±1.5mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	69N
Spinal Quincke needle - Size 19G, length 89 mm	
Size	19G (Ø1.10 mm)
Length	89 mm (±1.5mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	69N
Spinal Quincke needle - Size 20G, length 89 mm	
Size	20G (Ø0.90 mm)
Length	89 mm (±1.5mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	54N
Spinal Quincke needle - Size 21G, length 89 mm	
Size	21G (Ø0.80 mm)
Length	89 mm (±1.5mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	44N
Spinal Quincke needle - Size 22G, length 89 mm	
Size	22G (Ø0.70 mm)
Length	89 mm (±1.5mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	40N
Spinal Quincke needle - Size 23G, length 89 mm	
Size	23G (Ø0.60 mm)
Length	89 mm (±1.5mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	34N
Spinal Quincke needle - Size 24G, length 89 mm	
Size	24G (Ø0.55 mm)

ATTESTED

NOTARY PUBLIC
(INDIA)

24 FEB 2020
33675
21.05.2020

Length	89 mm (± 1.5 mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μ m
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	34N
Spinal Quincke needle - Size 25G, length 89 mm	
Size	25G ($\varnothing 0.50$ mm)
Length	89 mm (± 1.5 mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μ m
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	22N
Spinal Quincke needle - Size 26G, length 89 mm	
Size	26G ($\varnothing 0.45$ mm)
Length	89 mm (± 1.5 mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μ m
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	22N

Spinal Quincke needle - Size 27G, length 89 mm	
Size	27G ($\varnothing 0.40$ mm)
Length	89 mm (± 1.5 mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μ m
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	22N
Spinal Quincke needle - Size 28G, length 89 mm	
Size	28G ($\varnothing 0.36$ mm)
Length	89 mm (± 1.5 mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μ m
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	22N
Spinal Quincke needle - Size 29G, length 89 mm	
Size	29G ($\varnothing 0.33$ mm)
Length	89 mm (± 1.5 mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μ m
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	22N
Spinal Quincke needle - Size 30G, length 89 mm	
Size	30G ($\varnothing 0.30$ mm)
Length	89 mm (± 1.5 mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μ m
Blunting radius	no more than 0.03 mm

ATTESTED

NOTARY PUBLIC
(INDIA)

24 FEB 2020
21.03.2020

Needle head-tube connection stability	22N
Spinal Quincke needle - Size 22G, length 89 mm with introducer needle, size 18G, length 38 mm	
Needle size	22G (Ø0.70 mm)
Needle length	89 mm (±1.5mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	40N
Introducer needle size	18G (Ø1.20 mm)
Introducer needle length	38 mm (±1.5mm)
Spinal Quincke needle - Size 23G, length 89 mm with introducer needle, size 19G, length 38 mm	
Needle size	23G (Ø0.60 mm)
Needle length	89 mm (±1.5mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	34H
Introducer needle size	19G (Ø1.10 mm)
Introducer needle length	38 mm (±1.5mm)

Spinal Quincke needle - Size 24G, length 89 mm with introducer needle, size 20G, length 38 mm	
Needle size	24G (Ø0.55 mm)
Needle length	89 mm (±1.5mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	34N
Introducer needle size	20G (Ø0.90 mm)
Introducer needle length	38 mm (±1.5mm)
Spinal Quincke needle - Size 25G, length 89 mm with introducer needle, size 21G, length 38 mm	
Needle size	25G (Ø0.50 mm)
Needle length	89 mm (±1.5mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	22N
Introducer needle size	21G (Ø0.80 mm)
Introducer needle length	38 mm (±1.5mm)
Spinal Quincke needle - Size 26G, length 89 mm with introducer needle, size 21G, length 38 mm	
Needle size	26G (Ø0.45 mm)
Needle length	89 mm (±1.5mm)
Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	22N
Introducer needle size	21G (Ø0.80 mm)
Introducer needle length	38 mm (±1.5mm)
Spinal Quincke needle - Size 27G, length 89 mm with introducer needle, size 22G, length 38 mm	
Needle size	27G (Ø0.40 mm)
Needle length	89 mm (±1.5mm)

ATTESTED

NOTARY PUBLIC
(INDIA)

24 FEB 2020
21.05.2020

Bevel type	Quincke
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μ m
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	22N
Introducer needle size	22G ($\varnothing$ 0.70 mm)
Introducer needle length	38 mm ($\pm$ 1.5mm)
Name	Meaning
<i>Spinal Pencil Point needle - Size 16G, length 89 mm</i>	

Size	16G ($\varnothing$ 1.60 mm)
Length	89 mm ($\pm$ 1.5mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μ m
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	69N

Spinal Pencil Point needle - Size 17G, length 89 mm

Size	17G ($\varnothing$ 1.40 mm)
Length	89 mm ($\pm$ 1.5mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μ m
Blunting radius	no more than 0.03 mm

Needle head-tube connection stability	69N
---------------------------------------	-----

Spinal Pencil Point needle - Size 18G, length 89 mm

Size	18G ($\varnothing$ 1.20 mm)
Length	89 mm ($\pm$ 1.5mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μ m
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	69N

Spinal Pencil Point needle - Size 19G, length 89 mm

Size	19G ($\varnothing$ 1.10 mm)
Length	89 mm ($\pm$ 1.5mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μ m
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	69N

Spinal Pencil Point needle - Size 20G, length 89 mm

Size	20G ($\varnothing$ 0.90 mm)
Length	89 mm ($\pm$ 1.5mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μ m
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	54N

Spinal Pencil Point needle - Size 21G, length 89 mm

Size	21G ($\varnothing$ 0.80 mm)
Length	89 mm ($\pm$ 1.5mm)

ATTESTED

**NOTARY PUBLIC
(INDIA)**

33675
24 FEB 2020

Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	44N
Spinal Pencil Point needle – Size 22G, length 89 mm	
Size	22G ($\varnothing 0.70$ mm)
Length	89 mm (± 1.5 mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	40N


Spinal Pencil Point needle – Size 23G, length 89 mm	
Size	23G ($\varnothing 0.60$ mm)
Length	89 mm (± 1.5 mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	34N

Spinal Pencil Point needle – Size 24G, length 89 mm	
Size	24G ($\varnothing 0.55$ mm)
Length	89 mm (± 1.5 mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	34N

Spinal Pencil Point needle – Size 25G, length 89 mm	
Size	25G ($\varnothing 0.50$ mm)
Length	89 mm (± 1.5 mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	22N

Spinal Pencil Point needle – Size 26G, length 89 mm	
Size	26G ($\varnothing 0.45$ mm)
Length	89 mm (± 1.5 mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	22N

Spinal Pencil Point needle – Size 27G, length 89 mm	
Size	27G ($\varnothing 0.40$ mm)
Length	89 mm (± 1.5 mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 μm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	22N

ATTESTED

NOTARY PUBLIC
(INDIA)

24 FEB 2020

<i>Spinal Pencil Point needle – Size 28G, length 89 mm</i>	
Size	28G (Ø0.36 mm)
Length	89 mm (±1.5mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	22N


<i>Spinal Pencil Point needle – Size 29G, length 89 mm</i>	
Size	29G (Ø0.33 mm)
Length	89 mm (±1.5mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	22N

<i>Spinal Pencil Point needle – Size 30G, length 89 mm</i>	
Size	30G (Ø0.30 mm)
Length	89 mm (±1.5mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	22N

<i>Spinal Pencil Point needle – Size 22G, length 89 mm with introducer needle, size 18G, length 38 mm</i>	
Needle size	22G (Ø0.70 mm)
Needle length	89 mm (±1.5mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	40N
Introducer needle size	18G (Ø1.20 mm)
Introducer needle length	38 mm (±1.5mm)

<i>Spinal Pencil Point needle – Size 23G, length 89 mm with introducer needle, size 19G, length 38 mm</i>	
Needle size	23G (Ø0.60 mm)
Needle length	89 mm (±1.5mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	34N
Introducer needle size	19G (Ø1.10 mm)
Introducer needle length	38 mm (±1.5mm)

<i>Spinal Pencil Point needle – Size 24G, length 89 mm with introducer needle, size 20G, length 38 mm</i>	
Needle size	24G (Ø0.55 mm)
Needle length	89 mm (±1.5mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	34N
Introducer needle size	20G (Ø0.90 mm)
Introducer needle length	38 mm (±1.5mm)

ATTESTED

NOTARY PUBLIC
(INDIA)

24 FEB 2020


21.05.2020

Spinal Pencil Point needle – Size 25G, length 89 mm with introducer needle, size 21G, length 38 mm	
Needle size	25G (Ø0.50 mm)
Needle length	89 mm (±1.5mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	22N
Introducer needle size	21G (Ø0.80 mm)
Introducer needle length	38 mm (±1.5mm)
Spinal Pencil Point needle – Size 26G, length 89 mm with introducer needle, size 21G, length 38 mm	
Needle size	26G (Ø0.45 mm)
Needle length	89 mm (±1.5mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	22N
Introducer needle size	21G (Ø0.80 mm)
Introducer needle length	38 mm (±1.5mm)
Spinal Pencil Point needle – Size 27G, length 89 mm with introducer needle, size 22G, length 38 mm	
Needle size	27G (Ø0.40 mm)
Needle length	89 mm (±1.5mm)
Bevel type	Pencil Point
Rigidity	Rigidity 51-58 HRC
Roughness	Roughness no more than 0.16 µm
Blunting radius	no more than 0.03 mm
Needle head-tube connection stability	22N
Introducer needle size	22G (Ø0.70 mm)
Introducer needle length	38 mm (±1.5mm)

The needle tube - the introducer needle hub connection must withstand 10 N for needles with a nominal outer diameter of less than 0.6 mm and 20 N for needles with a nominal outer diameter of 0.6 mm or more.

Sizes and needle tip configuration are given in the technical documentation.
Table 2 – Needle color code

Nominal outer diameter of the needle, mm	color
0.30	Yellow
0.33	Red
0.36	Blue- Green
0.40	Light gray
0.45	Brown
0.50	Orange
0.55	Light violet

ATTESTED

NOTARY PUBLIC
(INDIA)


24 FEB 2020

33675
21.05.2020

0.60	Dark blue
0.70	Black
0.80	Dark green
0.90	Yellow
1.10	Cream
1.20	Pink
1.40	Red
1.60	White

Materials:

All materials of component parts comply with the specifications and do not contain any components of animal origin

No.	Material use site	Main material	Brand, manufacturer	Direct body contact (Yes / No)
1	Needle cap	Polypropylene	J 800-1, HY-OSUNG Corporation	No
2	Needle hub	Polypropylene	J 700-1, HY-OSUNG Corporation	No
3	Introducer needle hub	Polypropylene	J 700-1, HY-OSUNG Corporation	No
4	Needle / introducer needle	Stainless Steel	SS 304, Henan hzz iron and steel Co., Ltd	Yes
5	Dye (pigment) Yellow Red Blue Light gray Brown Orange Light violet Dark blue Black Dark green Yellow Cream Pink Red White	Manufacturer Focus Pigment Sdn. Bhd Yellow 1422G Red 5312P Blue 1513H Gray GR1 Brown 9640(P.O.64) Orange 9130 Violet 7230R (P.V.23) Blue 1535 (P.R.15 : 3) BlackV-RL Green GFP-H (P.G.7) Yellow HGR Yellow 8323H (P.Y.83) Red 1222 Red 5312P(P.R.53:1) -	 ATTESTED  NOTARY PUBLIC (INDIA)	No
6	Paper for individual packaging	Paper for medical purpose	OCPSxxxx10x, AMCOR Flexibles	No
7	Film for individual packaging	Polyethylene PP SM498 + Polypropylene PE HI	Polyethylene PP SM498 + Poly-	No

		11600	propylene PE HI 11600, AMCOR Flexibles	
8	Two-layer box	Carton	--	No
9	Shipping box	High strength wrapping paper	--	No
10	Ethylene oxide gas	30% EO + 70% CO ₂	--	No

10. PACKAGE

Spinal needles in various designs are available in individual sterile package that meets the requirements of ISO 11607-1, ISO 11607-2, and they are packed in a 2-layer secondary package. Clause 6 of GOST P ISO 11607 is mentioned in document No. TPM/BLSTR/0357.

The package is free of holes, cracks, tears, wrinkles, wrinkles, thickenings or thinning.

Package size 175x32 mm

Packing material density: 420 mm x 150 µm, paper – 410mm x 60(±2) g/m².

Tensile strength – not less than 12 MPa

Ultimate tensile strength: ≥ 66 N/15 mm,

Breathability – 200-400 cm³/m²•min.

Tear resistance – not less than 20 N.

Continuous structure of the adhesive layer.

Pore size: ≤ 20 µm

Tensile strength (dry): ≥ 250 kPa

Tensile strength ≥ 350 mN

Bonding density: 60 + 10 g/m²

Seal strength ≥ 90 cN / 15 mm

Bonding strength – not less than 0.12 MPa.

Adhesive layer width – not less than 8 mm.

Sterilization validity - 5 years.

11. STERILIZATION

Spinal needles in various designs are delivered sterile in a sealed package. They are sterilized by ethylene oxide. The sterilization process is validated in accordance with the requirements of ISO 11135:2014.

Spinal needles in various versions are for single use only. Re-sterilization and reuse are prohibited.

12. LABELLING

The consumer package shall contain the following information:

a) needle designation in accordance with:

- 1) nominal outer diameter of the needle tube in mm;
- 2) nominal length of needle tube in mm;

b) the word "STERILE" or the corresponding symbol;

c) batch code with the word "BATCH" or the corresponding symbol with the date of sterilization (month and year);

d) name and / or trademark of the manufacturer;

e) the words "valid until ..." (month and last two numbers of the year) or the corresponding symbol;

f) warning to check the integrity of consumer packaging before use or the corresponding symbol.

ATTESTED

NOTARY PUBLIC
(INDIA)

24 FEB 2020

21.05.2020

Group package

The group package (if any) shall contain the following information:

a) needle designation in accordance with Chapter 7 GOST ISO 7864-2011 and, if necessary, words: "thin-walled" or "super-thin-walled";

b) word "STERILE" or the corresponding symbol;

c) wording "FOR SINGLE USE ONLY" or the corresponding symbol.

Note - Additional information for applying a symbol is given in Annex C;

d) warning to check the integrity of consumer packaging before use or the corresponding symbol;

e) batch code or the corresponding symbol;

f) date (year and month) of sterilization.

Note - The first two numbers of the batch code may indicate the sterilization date;

g) the words "valid until ..." (month and last two numbers of the year) or the corresponding symbol;

h) name and address of the manufacturer or supplier;

i) storage and transportation conditions.

Shipping package

Shipping package shall contain the following information:

a) needle designation in accordance with

1) nominal outer diameter of the needle tube in mm;

2) nominal length of needle tube in mm, and, if necessary, words: "thin-walled" or "super-thin-walled";

b) batch code or the corresponding symbol;

c) word "STERILE" or the corresponding symbol;

d) date (year and month) of sterilization;

e) name and/or trademark of the manufacturer;

f) storage and transportation conditions.



Labelling symbols:


ATTESTED

NOTARY PUBLIC
(INDIA)

24 FEB 2020

	Do not re-use
	Caution! Consult accompanying documentation!
	Use before
	Batch code
	Catalog number
	Sterile, method of sterilization using ethylene oxide
	Do not re-sterilize
	Do not use if damaged package
	Authorized Representative in the European Community
	Non-pyrogenic
	Wording: Non toxic

ATTESTED

(Signature)
NOTARY PUBLIC
(INDIA)

The external cardboard box and shipping cardboard box are labeled using additional symbols on the proper transportation and storage conditions, such as:

24 FEB 2020

33075
21.05.2020

- Temperature limitation 
- Keep dry 
- Keep out of sunlight 
- Non pyrogenic 
- Non toxic

13. SAFETY REQUIREMENTS

- This medical device is intended to be used exclusively for its intended purpose and only by or under the supervision of a qualified licensed doctor or other medical practitioner as prescribed.
- Any medical methods and procedures described in certain publications do not cover ALL acceptable medical protocols; nor do they substitute the experience and judgment of a doctor in the treatment of each individual patient.
- Keep sterile.
- Do not use if package is opened or damaged.
- Do not use if there are any damage or defect.
- Observe general precautions to avoid exposure to bloodborne pathogens during use.

14. ENVIRONMENTAL REQUIREMENTS.

When used, transported and stored, spinal needles in various designs do not adversely affect humans and the environment. When using the product, special protection of personnel and patient is not required.

15. LIST OF INTERNATIONAL REGULATORY DOCUMENTS / STANDARDS THE MEDICAL DEVICE MEETS.

Reference	Name
EN ISO 13485:2016	Medical devices. Quality management systems. Requirements for regulatory purposes
BS EN 556-1:2001	Sterilization of medical devices. Requirements for medical devices to be designated "STERILE". Requirements for terminally sterilized medical devices
BS EN 868-5:2009	Packaging for terminally sterilized medical devices. Sealable pouches and reels of porous and plastic film construction. Requirements and test methods.
BS EN 980:2008	Medical devices. Symbols for use in the labelling of medical devices
BS EN 1041:2008	Information supplied by the manufacturer of medical devices
BS EN ISO 14971:2009	Application of risk management to medical devices

24 FEB 2020

ATTESTED

NOTARY PUBLIC
(INDIA)

Reference	Name
ISO 10993-1:2009/ Cor.1:2010	Biological evaluation of medical devices — Part 1: Evaluation and testing within a risk management process.
ISO 10993-2:2006	Biological evaluation of medical devices. Part 2: Animal welfare requirements
ISO 10993-4:2002/Amd.1:2006	Biological evaluation of medical devices. Part 4: Selection of tests for interactions with blood.
ISO 10993-7:2008/ Cor.1:2009	Medical devices. Biological evaluation of medical devices. Part 7: Ethylene oxide sterilization residuals
ISO 10993-10:2010	Biological evaluation of medical devices. Part 10: Tests for irritation and skin sensitization
ISO 10993-11:2006	Biological evaluation of medical devices. Part 11: Tests for systemic toxicity
ISO 10993-12:2012	Biological evaluation of medical devices. Part 12: Sample preparation and reference materials
ISO 10993-18:2005	Biological evaluation of medical devices. Part 18: Chemical characterization of materials
ISO 11135-1:2014	Sterilization of health care products. Ethylene oxide. Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices
EN ISO 11607-1:2019	Packaging for terminally sterilized medical devices. Part 1: Requirements for materials, sterile barrier systems and packaging systems
ISO 11607-2:2019	Packaging for terminally sterilized medical devices. Part 2: Validation requirements for forming, sealing and assembly processes
ISO 14971:2012	Medical devices. Application of risk management to medical devices. Amendment 1. Justification of requirements
ISO 15223-1:2016	Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied

24 FEB 2020

16. TRANSPORTATION, STORAGE AND OPERATION CONDITIONS

Spinal needles in various designs are transported by all means of transport in covered vehicles in accordance with the transportation requirements and applicable to vehicles of each type.

When transported in shipping package, they are resistant to climatic factors at temperature conditions from 0 to + 40°C, air humidity 26 - 65% and mechanical stress during vibration loads in the frequency range 10 - 55 Hz, movement amplitude 0.35 mm and shock acceleration with a peak acceleration of 10g with a duration of impact acceleration of 16 ms.

When stored in shipping package, they are resistant to climatic factors at temperature conditions from 0 to + 40°C, air humidity 26 - 65% (without condensation). Protect from direct sunlight, heat and freezing.

During operation, they are resistant to climatic factors for the temperature regime U6 from 32 to 42°C.

ATTESTED

(Signature)

NOTARY PUBLIC (INDIA) 33675 21.05.2020

17. WARRANTY

Poly Medicure Limited GUARANTEES THAT THIS PRODUCT IS MANUFACTURED IN ACCORDANCE WITH APPLICABLE STANDARDS AND SPECIFICATIONS. PATIENT'S CONDITION, CLINICAL TREATMENT AND SERVICE OF THE PRODUCT MAY AFFECT THE OPERATION OF THE PRODUCT. THIS PRODUCT SHOULD BE USED IN ACCORDANCE WITH THE RECOMMENDATIONS AND PRESCRIPTION OF A DOCTOR.

Guaranteed shelf life is 5 years from the production date, provided that the packaging is not damaged or opened.

18. SHELF LIFE

Sterility retention period is 5 years. After the sterility expiration spinal needles in various designs must be disposed of in accordance with the requirements of local authorities.

19. DISPOSAL

The products shall be disposed in accordance with the standards adopted in medical practice, as well as relevant local, state and federal laws and regulations, such as class B waste SanPin 2.1.7.2790-10 "Sanitary and epidemiological requirements for the management of medical waste".

18. Authorized representative in the Russian Federation:

Advamedics LLC

Room 15, Floor 2, building 2, 10, Korovinskoe Shosse, Moscow, 127486, Russia

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