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USER MANUAL

**TRIGEN SURESHOT electromagnetic targeting surgical system for
nails distal locking during trauma operations with accessories**

**Instructions for Care, Maintenance, Cleaning and Sterilization of
Smith & Nephew Orthopaedics Devices**

State of Tennessee)

County of Shelby)

Subscribed (or affirmed) before me this 25th day
of Feb. (month), 2021, by Maureen Whitson
(name of signer), whose identity was proven to me by personal
knowledge or satisfactory evidence.

Deborah H. Phillips

(Signature of Notary)

(Seal of Notary)

Total number of pages:

Name: Maureen whitson Mall

Title: Regulatory Affairs Specialist

TRIGEN SURESHOT distal targeting system

Operations/Service Manual

TRIGEN SURESHOT Distal Targeting System V4.0


**TRIGEN
SURESHOT**
Distal Targeting System

CECE 2797

81103554-A

TRIGEN SURESHOT Trauma Interface and Targeter (2nd Generation)

TRIGEN SURESHOT Trauma Interface

Glossary of Symbols


Symbol	Value	Symbol	Value
	Medical Device		Humidity limit
	Electronic instructions for use		Atmospheric pressure
	Power ON (connection to the mains)		Temperature range
	Power OFF (disconnection from the mains)		Follow instructions for use
	Sterilized using ethylene oxide		Use-by date
	Sterilized using irradiation		Batch number
	Equipment classification – Patient Isolation/Type BF Applied Part		General Warning Sign
	Fuse		Date of manufacture
	CE Mark		Unique Device Identifier
	Authorized representative in the European community		Serial number
	Testing Lab certification		Catalog number
	Manufacturer		EU Not for general waste
	Unsterile		Do not reuse
	Prescription use only		Conformity mark for mandatory certification of products in the Russian Federation
	Consult instructions for use		

Table of contents

Device Description	7
How Supplied	7
Intended Users	7
Clinical Benefits	7
Intended Use	7
Indications for Use	7
Contraindications	7
Target Patient Population	7
Training	7
Plausibility Check	7
Warnings	8
Precautions	8
Accessibility of Documentation	8
The Operator Checks and Decides	8
Responsibility of Smith & Nephew	8
Cleaning and Sterilization of Instruments	8
Repair or Modifications to the System	8
Modifications/Additions to the Software	8
IR Safety	9
Maintenance	9
Recycling	9
Possible Adverse Events	9
Adverse Event Reporting	9
Device Lifetime	9
Quick Start Guide	10
Where to Find the Components	10
Front Panel Layout	10
Rear Panel Layout	10
Applied Parts for Targeting	11
Surgical Procedure – O.R. Preparation	12
Trauma Interface Setup	12
Software Configuration	12
TRIGEN SURESHOT Targeter Connection	13
Surgical Procedure – After IM Nail Assembly to the Drill Guide	15
Probe Selection and Assembly	15
Probe Connection	16
Drill Sleeve Selection	16
Implant Selection Screen	17
Drill Sleeve Attachment	17

Locking Hole Accuracy Check in the Operative Field	18
Optional – Field Accuracy Check	18
On-Screen Nail Rotation	18
Rotation Using Targeter	19
Rotation Using Monitor	19
Surgical Procedure – Distal Targeting	19
Probe and Set Stop Insertion	20
Skin Incision	20
Targeting the Locking Hole	20
Surgical Procedure – Drill Depth Measurement	21
Option 1: Measure off Drill Bit Calibrations	21
Option 2: Measure Using Drill Depth Software	21
Screw Insertion	23
Using the Combo Drill Sleeve	23
Navigation Screen Operation	23
Overview Mode	23
Drilling Mode	23
Shutting Down the System, Storage, and Transport	23
Shutting Down the System	24
Storage and Transport	24
Maintenance	24
Field Accuracy Check	24
Field Accuracy Check Steps	24
Software Upgrades	25
Cleaning of the TRIGEN SURESHOT Trauma Interface	25
.....	25
Replacing Fuses	25
Technical Information	26
Service	26
TRIGEN SURESHOT Distal Targeting System	26
Electrical Interference	26
Environmental Protection	26
Troubleshooting	27
Technical Specifications	30
Electrical Safety	30
Parts and Accessories List	30
Guidance and Manufacturer's Declaration	32
Electromagnetic Emissions, Immunity and Separation Distances	32
Compliant Cables and Accessories	32
Safe Disposal	36
For Further Information	36

The following names used throughout this document are identical to the name of the medical device **TRIGENSURESHOT electromagnetic targeting surgical system components for nails distal locking during trauma operations with accessories:**

- System;
- TRIGEN SURESHOT System;
- TRIGEN SURESHOT Distal Targeting System;
- Instruments;
- Devices.

Device Description

The TRIGEN SURESHOT Distal Targeting System is a software-controlled electromagnetic tracking system. It assists the surgeon in locating and positioning screws in an intramedullary nail implant during orthopedic trauma surgery.

How Supplied

The TRIGEN SURESHOT Distal Targeting System devices are packaged separately, with the exception of some reusable surgical devices that may be provided as a set within an instrument tray.

Symbols on the package labeling are used to indicate if each device is provided cleaned and sterile or non-sterile. Devices provided sterile are sterilized to a Sterility Assurance Level (SAL) of 10^{-6} . Sterile packaged devices are supplied in protective sterile barrier packaging. Inspect packages for punctures or other damage prior to surgery. If the sterile barrier has been broken, return the device to Smith & Nephew, Inc.

With the exception of the TRIGEN SURESHOT Trauma Interface, devices provided non-sterile must be cleaned and sterilized before use according to the information provided in the Smith & Nephew document *Instructions for Care, Maintenance, Cleaning and Sterilization of Smith & Nephew Orthopedics Devices*, available on the Smith & Nephew website. The TRIGEN SURESHOT Trauma Interface should be cleaned before use, per the *Cleaning of the TRIGEN SURESHOT Trauma Interface* section.

Intended Users

Smith & Nephew's TRIGEN SURESHOT Distal Targeting System is intended for use by physicians who have thorough knowledge of the anatomy, physiology, and surgical principles involved.

Physicians are strongly encouraged to obtain instruction from experienced clinicians or to observe surgical application of the system prior to initial use.

Clinical Benefits

The TRIGEN SURESHOT system is expected to have the following clinical benefits:

- Shortened procedure time.
- Reduced radiation exposure.

Intended Use

The TRIGEN SURESHOT Distal Targeting System is intended to be an intraoperative image guided localization system. It is a computer assisted orthopedic surgery tool to aid the surgeon with drill positioning for screws during intramedullary nail implantation. It provides information to the surgeon that is used to place surgical instruments during surgery utilizing intraoperatively obtained electromagnetic tracking data.

Indications for Use

The TRIGEN SURESHOT Distal Targeting System is indicated for long bone fractures treated with intramedullary nails in which the use of stereotactic surgery may be appropriate.

Contraindications

The screw targeting software application for the TRIGEN SURESHOT System is contraindicated for all IM nails other than Smith & Nephew TRIGEN META-NAIL, TRIGEN META-TAN, TAN, FAN Nails, Humeral, Pediatric and Adolescent nails. Do not operate the TRIGEN SURESHOT Targeter within 200 mm of an installed pacemaker. The magnetic field produced by the Targeter may interfere with the operation of the pacemaker.

Target Patient Population

No specific patient population is defined for TRIGEN SURESHOT Distal Targeting System. Rather, the target patient population for the system is aligned to the intended patient population for the corresponding intramedullary (IM) nails manufactured by Smith & Nephew.

Training

Only trained operators are allowed to use the TRIGEN SURESHOT System. The various operating instructions must be fully read and understood as part of the training. If any part of the instructions is not clear, please contact your local representative.

Plausibility Check

As with all technical equipment, malfunctions may occur due to improper use or, more rarely, technical failure. To reduce the risks involved with such technical malfunction the operation can be completed using manually controlled instruments, providing the malfunction is detected without delay.

It is, therefore, important to check the plausibility of the steps, as indicated by the system, and to carry out verification of the software targeting, particularly when using the system for the first time. Should there be any doubt regarding correct functioning, the targeting should be verified or a switch made to a traditional X-Ray technique.

Warnings:

- This device is provided non-sterile and must be cleaned and sterilized per *Instructions for Care, Maintenance, Cleaning and Sterilization of Smith & Nephew Orthopaedic Devices*.
- It is the health care professional's responsibility to be familiar with the appropriate surgical techniques prior to use of this device .
- Read these instructions completely prior to use .
- The patient should be warned of surgical risks, and made aware of possible adverse effects . It is recommended that the operating surgeon record all information given to the patient .
- The maximum temperature of the Targeter body can reach 47° Celsius at ambient room temperatures above 30° C . The Targeter body should not remain in constant contact with a patient's exposed skin for more than one minute .
- Verify there are no other metal objects (including metal triangles) in the immediate targeting area . Metal interference will cause the system to be inaccurate .
- To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth (=ground) .

Avoid spilling water or other liquids on electronic/electrical equipment . Use only Smith & Nephew disposables and accessories with the TRIGEN SURESHOT Distal Targeting System .

- Do not operate the TRIGEN SURESHOT Distal Targeting System within 200 mm of an installed pacemaker .
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Smith & Nephew TRIGEN SURESHOT Distal Targeting System, including cables specified by the manufacturer . Otherwise, degradation of the performance of this equipment could result ."
- Verify there are no other metal objects (including metal triangles) in the immediate targeting area . Metal interference will cause the system to be inaccurate .
- Do Not Resterilize SURESHOT Probes .
- The standard TRIGEN SURESHOT Drill Bits are made from magnetic stainless steel that will cause interference with the system and cannot be used . Be sure to use the gold non-magnetic TRIGEN SURESHOT Drill Bit or long golden non-magnetic drill bit .
- **STERILE DEVICES ONLY:** Contents are sterile unless the package is opened or damaged. **DO NOT RESTERILIZE.** For single use only. Discard any open unused product. Do not use after the expiration date .

Precautions

U.S. Federal law restricts this device to sale by or on the order of a physician.

Prior to each use, inspect the device to ensure it is functioning properly and not damaged. Do not use a damaged device.

Accessibility of Documentation

Please ensure that all instructions are kept in an easily accessible place for operating personnel.

The Operator Checks and Decides

All the information provided by the TRIGEN SURESHOT Distal Targeting System is to help the operator make decisions during the operation. The operator must check all the suggestions made by the system and is responsible for all decisions taken.

Responsibility of Smith & Nephew

In the event of improper use, Smith & Nephew accepts no responsibility or liability whatsoever for the functioning or utility of the TRIGEN SURESHOT Distal Targeting System when used in the operating theater.

Cleaning and Sterilization of Instruments

All instruments must be cleaned and sterilized before use. Detailed information on the cleaning and sterilization of reusable components is contained in the separate Instructions for Care, Maintenance, Cleaning and Sterilization of Smith & Nephew Orthopedics Devices available on the Smith & Nephew website.

Repair or Modifications to the System

The user is not permitted to modify or service the equipment. There are no serviceable parts inside the unit. Refer all service to authorized personnel.

Modifications/Additions to the Software

The user is not permitted to install or uninstall software. Any new software must be installed by the manufacturer or by authorized personnel.

It is only allowed to connect equipment to the interface and power supply connections of the TRIGEN SURESHOT Distal Targeting System which are IEC 60601-1 approved and which are approved by Smith & Nephew Orthopaedics. Do not modify this equipment without authorization of the manufacturer.

MR Safety

The TRIGEN SURESHOT Distal Targeting System does not include implantable devices and has not been evaluated for use in the Magnetic Resonance (MR) environment. Consult the respective instructions for use of the compatible intramedullary nails for information regarding the compatibility of the implantable devices with MR environments.

Maintenance

To verify accurate functionality, the device should be checked per the Maintenance section of this document. This accuracy check must be performed at least once every 12 months. If this accuracy check is not performed as defined in the previous paragraph, all warranty claims expire and the device is operated at the user's own risk.

Recycling

Old electrical and electronic equipment must be disposed of separately and may not be included in the regular domestic waste. Alternatively, the unit can be handed over to Smith & Nephew Orthopaedics for correct recycling.

Possible Adverse Events

Users should observe the following possible adverse events:

1. Loosening, bending, cracking or fracture of the implant components or instruments
2. Infections, both deep and superficial.
3. Tissue reactions which include macrophage and foreign body reactions adjacent to implants
4. Irrational injury of soft tissues, including impingement syndrome.
5. Screw disengagement leading to increase surgery time, implant instability, and/or revision surgery.
6. Increased surgery time due to missing the targeted hole, software malfunction/errors, device incompatibility or non-recognition, instrument breakage, or instrument malfunction.

Additional adverse events associated with use of the compatible intramedullary nails are available in the respective instructions for use.

Adverse Event Reporting

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer (complaints@smith-nephew.com) and to your national authority.

Device Lifetime

All hardware (except disposables components) are designed for a lifetime of at least 10 years. The expiration date of drills and probes is 10 years.

The TRIGEN SURESHOT Distal Targeting System includes single-use devices and reusable devices. The lifetime of single-use devices is the single use of the device. The TRIGEN SURESHOT Trauma Interface end of life is determined by the end of the serviceable life. See the Plausibility Check, Maintenance and Service sections of this document for additional information regarding the functional checks and serviceability of the trauma interface.

The useful life of the TRIGEN SURESHOT Distal Targeting System reusable surgical instruments depends on many factors including the method and duration of use, and the handling between uses. Careful inspection and function test of the device before use is the best method of determining the end of serviceable life for the medical device. Do not use devices that show evidence of damage and wear. Evidence of damage and wear on a device may include but is not limited to corrosion (i.e. rust, pitting), discoloration, excessive scratches, flaking, wear and cracks. Improperly functioning devices, devices with unrecognizable marking, missing or removed (buffed off) part numbers, damaged and excessively worn devices should not be used.

⚠ WARNING: Do not unplug the power while the system is running!

⚠ WARNING: Danger of damage and tipping over!

Note: Place the unit on a firm, level surface capable of holding at least 10 kg (22 lbs).

⚠ WARNING: To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

Quick Start Guide

Where to Find the Components

1. Trauma Interface
2. Power cord (country-specific)
3. User Manual slot

Front Panel Layout

1. Probe sensor ports
2. Targeter port
3. Touch screen interface

Rear Panel Layout

1. Power switch
2. Power connection (Appliance Inlet)
3. USB connections (service only)
4. VGA video output (1280 x 800 resolution)
6. HDMI video output (1280 x 800 resolution)

Note: External equipment intended for connection to signal input, signal output or other connectors shall comply with the relevant product standard e.g. IEC 60950-1 or IEC 62368-1 for IT equipment and the IEC 60601-series for medical electrical equipment. In addition, all such combinations – Medical Electrical Systems – shall comply with the safety requirements stated in the collateral standard IEC 60601-1-1 or the general standard IEC 60601-1, edition 3, clause 16. Any equipment not complying with the leakage current requirements in IEC 60601-1 shall be kept outside the patient environment i.e. at least 1.5 m from the patient support or shall be supplied via a separation transformer to reduce the leakage currents.

Any person who connects external equipment to signal input, signal output or other connectors has formed a Medical Electrical System and is therefore responsible for the system to comply with the requirements. If in doubt, contact a qualified medical technician or your local representative.

Applied Parts for Targeting

TRIGEN SURESHOT Targeter
(Body Applied Part Type BF)

WARNINGS:

- The maximum temperature of the Targeter body can reach 47° Celsius at ambient room temperatures above 30° C . The Targeter body should not remain in constant contact with a patient's exposed skin for more than one minute .
- Verify that the Targeter housing is not damaged (holes, tears, cracks) . If the housing or the connector is damaged, the Targeter is no longer safe to use .
- This device is provided non-sterile and must be cleaned and sterilized per Instructions for Care, Maintenance, Cleaning and Sterilization of Smith & Nephew Orthopaedic Devices prior to use.

Cautions

- The Targeter will be operated within the sterile field and may have contact with the skin of the patient. The drill sleeve inserts will be used in the incision and have direct bone contact.
- If the Targeter is not recognized after its connection to the system, the Targeter is defective and must be exchanged. (See also instrument connection.)
- Broken or damaged instruments must be exchanged immediately and sent back to Smith & Nephew, Inc.

TRIGEN SURESHOT Probes (Body Applied Part Type BF)

TRIGEN SURESHOT META-NAIL Standard Drill Guide Probe
The probe will be used as an intramedullary tool inside the nail placed in the patient's bone.

WARNINGS:

- This device is provided sterile by ethylene oxide gas and is single use
- In case the packaging is compromised, please replace with new sterile packaged probe
- Do Not Resterilize SURESHOT Probes .

CAUTION: If the probe is not recognized after its connection to the system, the probe is defective and must be exchanged. (See also instrument connection).

CAUTION: Broken or damaged instruments must be exchanged immediately and sent back to Smith & Nephew, Inc.


Surgical Procedure – O .R .

Preparation

This procedure will cover only the specific steps of freehand targeting of intramedullary locking holes using the TRIGEN SURESHOT Distal Targeting System. For the full surgical procedure, please refer to the specific surgical technique for the TRIGEN IM Nail System being implanted.

Trauma Interface Setup

After the sterile areas have been established, place the Trauma Interface in the desired non-sterile location and turn on the power switch.

CAUTION: No other electrical devices should be placed near the Trauma Interface. See the “Guidance and Manufacturer’s Declaration – Separation Distances” table at the end of this document.


Note: If the Trauma Interface does not power on, make sure the switch is in the “on” position.

Note: The means for mains disconnection of power to the Trauma Interface is the appliance inlet located below the power switch. The Trauma Interface should not be positioned such that it is difficult to reach this location.

Software Configuration

To enable or disable specific software features, use the following steps to make the desired changes. These software features could include, but are not limited to, TAN/FAN nails, the Japan nail option, META-TAN nails, the Drill Depth Measurement option, on screen nail rotation via the targeter, the Combo Drill Sleeve option, and shut down V2.


1. Press the “Menu” button.
2. Select the “About” button from the dropdown.


3. Select the “Administration” button.


- Press the "Settings" button in the upper right corner of the screen.

- Select "Language" on the left to bring up the following screen:

- Select the desired language from the box on the right.
- Select "Application" on the left to bring up the following screen:


- Press the feature to be changed to enable or disable that feature. Enabled features appear in black text with a filled check box; disabled features appear in gray text with an unfilled check box.
- Press the "Return" button in the upper left corner when finished.

- TRIGEN SURESHOT Distal Targeting System button in the center of the screen to return to the targeting application.


TRIGEN SURESHOT Targeter Connection

When the display prompts for tool connections, connect the TRIGEN SURESHOT Targeter (7169-2801) to the Targeter port on the Trauma Interface.

The SURESHOT-Targeter will change color to orange upon successful connection to the Trauma Interface.

CAUTION: The Targeter body may have contact with the patient and must remain in the sterile field at all times. Only the cable and connector may be removed from the sterile field.

CAUTION: Connect the Targeter at least 10 minutes prior to targeting in order to ensure proper accuracy.

Notes:

- When oriented as shown, the connector should assemble easily . Do not force the connector into the port .
- If the Targeter is properly connected to the system and the application remains in this screen for more than 30 seconds, the Targeter may have been damaged during cleaning/sterilization . In this case another Targeter has to be used .
- It is possible at any time to disconnect and reconnect tools when the application is running . The display will show a screen reporting the missing instrument .


Surgical Procedure – After IM Nail Assembly to the Drill Guide

Probe Selection and Assembly

Assemble the appropriate probe and set the stop for the TRIGEN IM Nail and Drill Guide that will be used.


CAUTION: Proper orientation of the probes to the set stops as shown is required. Failure to do so may result in inaccurate targeting.

TRIGEN SURESHOT
META-NAIL
Standard
Drill Guide Probe
(7169-2814)

Red probe
Use only with
META-NAIL
Standard Drill Guide
(7165-4502)

TRIGEN SURESHOT
META-NAIL
Semi-extended
Drill Guide Probe
(7169-2803)

Blue probe
Use only with
META-NAIL
Semi-extended Drill Guide
(7165-4524)

TRIGEN SURESHOT
Percutaneous
TAN/FAN
Drill Guide Probe
(7169-2815)

Green probe
Use only with Percutaneous
TAN/FAN nails , Drill Guide
(7163-1021)
and the META-TAN Drill Guide
(7163-1451)

TRIGEN SURESHOT
Humeral Nail
Drill Guide
Probe
(7169-1152)

Dark Green probe
Use only with
Humeral Drill
Guide (7175-1129)

Set the probe to the length of the TRIGEN IM Nail.

For META-NAIL,
notches should face anteriorly

For TAN/FAN nails and META-TAN
notches should face medially

For Humeral,
notches should face laterally

Probe Connection

Connect the probe to either of the probe sensor ports on the Trauma Interface.

Confirm that the tool connection has been verified.

Notes:

- When oriented as shown, the connector should assemble easily. Do not force the connector into the port. Simply try rotating the connector until the keys are oriented in the 12 o'clock position.
- If the probe is properly connected to the system and the application reports "Probe not found" for more than 10 seconds, the probe may be damaged or defective. In this case, the probe has to be exchanged.
- It is possible at any time to disconnect and reconnect tools when the application is running. The display will show a screen reporting the missing instrument.


When a humeral probe is connected, the following warning message will appear:

WARNING: Do not operate the TRIGEN SURESHOT Targeter within 200 mm of an installed pacemaker.


Drill Sleeve Selection

Select the length of the drill sleeve (7169-2804, 7169-2805, or 7169-1165 and 7169-1166) that will be used.

A different sleeve can be selected at any time during the procedure by choosing the drill sleeve option from the drop down menu OR by pressing the drill sleeve tab located in the lower right corner of the navigation screen.


Note: For humeral nails, only one drill sleeve option will appear (7169-1154).




Implant Selection Screen

Select the TRIGEN IM Nail and size that will be used.



META-NAIL


TAN/FAN

META-TAN

Humeral Nail

Note: A different TRIGEN IM Nail and/or size can be selected at any time during the procedure by choosing the Implant option from the drop down menu OR by pressing the implant tab located in the lower left corner of the navigation screen.

Drill Sleeve Attachment

Tightly secure the selected drill sleeve to the Targeter.

Note: The short, long, and long combo drill sleeves (7169-


2805, 7169-2804, and 7169-1166) can be loosened from the Targeter using the slot in the TRIGEN Hammer (7167-4082).

Note: The Humeral Drill Sleeve (7169-1154) can be loosened from the targeter using the slot in the TRIGEN Slotted Hammer (7175-1135).


Locking Hole Accuracy Check in the Operative Field

Insert the probe with the assembled set stop through the drill guide and cannulation of the TRIGEN IM Nail.

Attach the set stop to the drill guide, ensure that the probe is oriented correctly and that the set stop position and IM nail length match.

Guide (7165-4503) or the Standard TRIGEN™ Drill Guide (7163-1134).

CAUTION: All tool cables should be uncoiled completely and any excess cables should be kept out of the Targeter measurement volume.

CAUTION: To guarantee system accuracy, the accuracy check has to be performed directly in the operative field.

Place the IM nail with inserted probe directly beside the patient's limb. Position the Targeter over the locking holes and verify the accuracy on the Trauma Interface.

Optional – Field Accuracy Check

An optional field accuracy check can be performed at this point using the instructions in the Maintenance section.

For META-NAIL, notches should face anteriorly

Use only with the META Standard (7165-4502) or Semi-extended Drill Guide (7165-4524)

For TAN/FAN and META-TAN notches should face medially

Use only with the Percutaneous TAN/FAN Drill Guide (7163-1021) and the META-TAN Drill Guide (7163-1451)

For Humeral, notches should face laterally

Use only with the Humeral Drill Guide (7175-1129)


CAUTION: The TRIGEN™ SURESHOT™ Distal Targeting System cannot be used with the TRIGEN META-NAIL Extension Drill

On-Screen Nail Rotation


Rotation Using Targeter

One of the new features to version 4.0 software is the surgeon's ability to manipulate nail rotation by rotating the Targeter in hand as opposed to having to touch the Interface monitor.

1. To activate the on-screen nail rotation feature, pull Targeter away from the implant.

2. Once the Targeter is pulled back, the on-screen message will appear, instructing the operator to rotate Targeter to trigger nail rotation.

3. Now that the nail rotation feature is activated, the operator can rotate the Targeter to position the implant to appropriate patient orientation.


Note: The Targeter should remain constant. Moving the Targeter toward the implant while performing implant rotation will deactivate the nail rotation feature.

Note: This feature can be disabled in the Settings menu of the Software as described in the software configuration section of this manual.

Rotation Using Monitor

Each IM nail has several predefined views that are automatically selected depending on the position of the Targeter to the IM nail. Depending on the operating environment, these predefined views might not be appropriate and can be manually adjusted.

To rotate the view

Touch the screen near the outside and “drag” the view in a clockwise or counterclockwise direction.

To reset the view

The default view settings can be restored by touching the center.

Surgical Procedure – Distal Targeting

WARNING: Verify there are no other metal objects (including metal triangles) in the immediate targeting area. Metal interference will cause the system to be inaccurate.


Probe and Set Stop Insertion

Insert the probe and set stop into the IM nail as described in the locking hole accuracy check on the previous page.

Note: Reference the cautions in the “Locking hole accuracy check” section.

Skin Incision

Position the serrated tip of the drill sleeve (represented by the green circle) over the desired locking hole and make the initial incision.


Note: If performing this procedure using any of the iRIGEN femoral antegrade nails (TAN/FAN /Pediatric/Adolescent), the Anteversion Locking Guide (7169-2816) should be used along with a 4.0 mm Drill Sleeve/Drill Bit to keep the IM nail from rotating.


Targeting the Locking Hole

With the appropriate length TRIGEN SURESHOT 4.0 mm Drill Bit (7169-2810 or 7169-2811) or the TRIGEN SURESHOT Humeral AO Drill Bit (7169-1155) inserted into the Targeter, insert the tip of the drill sleeve (represented by the green circle) through the incision and down to bone.


Adjust the trajectory of the Targeter (represented by the red circle) until both the green and red circles are concentric and drill through to the far cortex.

WARNING: The standard SURESHOT Drill Bits are made from magnetic stainless steel that will cause interference with the system and cannot be used. Be sure to use the golden non-magnetic TRIGEN SURESHOT Drill Bit or long golden non-magnetic drill bit.

Note: The green ring must be fully within the hole of the IM nail displayed on the navigation screen to ensure accurate drilling.

When the Targeter is out of the preferred range or there is metal or electrical interference, the green and red Targeter circles on the navigation screen may become unstable and/or a warning message will be displayed. If the interference is excessive, the IM nail image on the navigation screen will disappear.

A Field Quality meter is located at the bottom of the screen. There are three field quality states: OK, Warning, and Error:

	<p>OK: Field Quality is adequate for proper function.</p>
	<p>Warning: Possible measurement error due to metal or electrical interference.</p>
	<p>Error: Interference is excessive. IM nail tracking not possible.</p>

If interference cannot be avoided, a standard X-Ray technique must be used.

CAUTION: All tool cables should be uncoiled completely and any excess cables should be kept out of the Targeter measurement volume.

Surgical Procedure – Drill Depth Measurement

Option 1: Measure off Drill Bit Calibrations

Begin drilling the hole at the near cortex. Before drilling through the far cortex, measure the length using the calibrations on the TRIGEN SURESHOT 4.0 mm Drill Bit (7169-2810 or 7169-2811), or the TRIGEN SURESHOT Humeral AO Drill Bit (7169-1155).


After measuring the length on the drill bit prior to drilling the far cortex, add the appropriate screw length in order to account for the thickness of the far cortex. Ensure the Targeter and the sleeve are pushed against the bone.

Option 2: Measure Using Drill Depth Software


When software is enabled for performing drill depth measurements, the following message will be displayed at the top of the screen, "Pull Targeter back fully to measure depth."

Begin drilling the hole at the near cortex. Before drilling through the far cortex, pull the Targeter back to the drill connection while keeping the drill bit in place. This action will trigger the software to display the Drill Depth Measurement window.



Note: Do not excessively rotate the Targeter while sliding it along the drill bit as this movement is used to cancel the Drill Depth Measurement software feature.

If the Targeter distance pulled back does not match the distance pushed back to bone then a warning message is displayed.


Push the Targeter drill sleeve forward to the bone to verify the length measured. This action verifies the drill bit did not change position during the Targeter movement. If the distance the Targeter was pulled back equals the distance pushed back to bone a confirmation message indicating an accurate measurement is displayed.

Note: After measuring the length on the drill bit prior to drilling the far cortex, add the appropriate screw length in order to account for the thickness of the far cortex.

Note: The Drill Depth Measurement software feature can be cancelled at any time by rotating the Targeter about the drill bit axis.

After a successful Drill Depth Measurement is performed the software will display the measured length at the drilled location.


Note: A Drill Depth Measurement software operation can be repeated on a drilled hole at any time. The newest measurement will be displayed at that location.

The user can delete a measurement at any time by tapping the measurement number on the touchscreen. A button will appear allowing the user to delete the measurement.

Screw Insertion

Using the TRIGEN SURESHOT Hexdriver (7169-2809 or 7169-1153), the screw may be inserted using the Targeter.

WARNING: The standard TRIGEN Hexdrivers are made from magnetic stainless steel that will cause interference with the system and cannot be used.

Using the Combo Drill Sleeve

TRIGEN SURESHOT Inner Drill Sleeve-Long (7169-1165)

TRIGEN SURESHOT Outer Screw Sleeve-Long (7169-1166)

1. While holding the Targeter to bone, remove the SURESHOT Inner Drill Sleeve-Long (7169-1165) from the SURESHOT Outer Screw Sleeve-Long (7169-1166).
2. Use the TRIGEN SURESHOT Hexdriver (7169-2809 or 7169-1153) to insert the appropriate length screw through the SURESHOT Outer Screw Sleeve-Long.
3. When screw insertion is complete, remove the probe from the IM nail.


Navigation Screen Operation

Overview Mode

When the Targeter is greater than 5 cm from the interlocking holes, the navigation screen will display the IM nail in the overview mode. This provides the user with a larger field of view in order to help find the general location of the interlocking holes.

CAUTION: If the desired hole in the nail is not oriented to be viewed on the screen, then the orientation of the Targeter is not aligned closely enough to that specific hole. Try holding


the Targeter generally anterior-posterior (AP) or generally mediolateral (ML) to orient the desired view.

Drilling Mode

When the Targeter is moved within 5 cm of the interlocking holes, the navigation screen will display the IM nail in the drilling mode. This provides the user with a smaller field of view that automatically zooms in to the interlocking holes.

The white lines displayed on either side of the IM nail can be used for targeting blocking screws. These lines are located 2.5 mm from the side of the IM nail for all IM nails 10 mm or larger in diameter. These lines are located 2 mm from the sides of 8.5 mm IM nails.

Shutting Down the System, Storage, and Transport

Shutting Down the System

1. Shut down the Trauma Interface by selecting the "Exit" option from the on-screen "Menu" button. Wait for the system to display the following screen:

2. Switch off the main power switch on the rear of the unit.
3. Unplug the power supply cable.

CAUTION: The system **MUST** fully shutdown (as shown above) before switching the power off.

CAUTION: The internal fan of the Trauma Interface runs continuously whenever the rear power switch is in the ON position. Switch power off whenever the unit is not in use.

Storage and Transport

- Coil the power cord and store in the transportation case.
- Consider the temperature and humidity range for storage and transportation on the transportation case.
- Place the Trauma Interface 2nd Generation (screen facing inwards) in the transportation case.

CAUTION: Do not place any objects on top of the screen or system when the Trauma Interface is stored in the shipping container. Damage to the screen or system may occur.

Maintenance

Field Accuracy Check

A field accuracy check procedure should be performed at least once a year or whenever the accuracy of a TRIGEN SURESHOT Probe or TRIGEN SURESHOT Targeter needs to be verified. This procedure can also be performed during surgery to verify all components are working correctly prior to their use on a patient.

CAUTION: This step should be performed at least once a year to ensure that the device is working properly.

Field Accuracy Check Steps



1. Attach TRIGEN SURESHOT Field Accuracy Gauge (7169-2808) to the TRIGEN SURESHOT Targeter. The knob on the Field Accuracy Gauge should be hand-tightened only.

2. Attach the META Set Stop (7169-2806) to the end of the Field Accuracy Gauge, insert a TRIGEN SURESHOT probe into the set stop, and set the depth to the "REF" mark on the probe body.


3. From the software "Menu" button located in the upper left corner of the screen, select the "Field Check" option.

4. A software window will appear informing the user if the TRIGEN SURESHOT Targeter and Probe combination is within the predefined accuracy parameters ("Pass" or "Fail" message).

If the field accuracy check fails, check the "Troubleshooting" section of this document for possible solutions.

Software Upgrades

Software upgrades are performed using a Smith & Nephew, Inc.- supplied USB memory stick and a Work Instruction document. The Work Instruction document should be filled out and returned to Smith & Nephew Customer Service upon completion of the upgrade.

Cleaning of the TRIGEN SURESHOT Trauma Interface

The TRIGEN SURESHOT Trauma Interface is used in the non-sterile area of the operating theater and is cleaned by wiping the surface with a clean, lint-free cloth and using detergent and water. The detergent should be a low foaming detergent with a pH between 6.0 and 8.0.

CAUTION: Care must be taken not to allow any liquid to pass into any electrical connections or the interior of the unit. Let the surfaces dry thoroughly before plugging in the unit. DO NOT steam sterilize the system. DO NOT submerge the system for any reason.

All other reusable instruments shall be cleaned and sterilized according to the *Instructions for Care, Maintenance, Cleaning*

and Sterilization of Smith & Nephew Orthopaedic Devices.

WARNINGS:

- Observe the manufacturer's instructions for dilution, exposure time, etc .
- DO NOT sterilize with ethylene oxide gas or steam sterilize the Trauma Interface .
- Probes are NOT reusable .

Replacing Fuses

See the "Technical Specifications" section for replacement fuse type.


1. Disconnect main power cord.
2. Use a screwdriver to open the fuse compartment door on the AC receptacle.

3. Pull out fuse carriers to exchange fuses.

4. Reinsert fuse carriers using the arrows on the inside of the fuse compartment door as a guide.

5. Snap the fuse compartment door closed.

WARNING: To avoid fire hazard, use only fuses of the correct type, voltage rating, and current rating .

CAUTION: Disconnect the power cord before exchanging any fuse.

Technical Information

Smith & Nephew will make available on request a list of all repairable exterior parts with descriptions. Interior schematics and circuit diagrams will be made available to qualified personnel only.

CAUTION: Service can only be performed by authorized Smith & Nephew personnel.

Service

TRIGEN SURESHOT Distal Targeting System

There are no user-serviceable components inside the TRIGEN SURESHOT Trauma Interface. Repairs and adjustments are to be performed only by Smith & Nephew authorized service centers.

If service becomes necessary, call an authorized Smith & Nephew customer service representative prior to returning the device and request a Return Authorization (RA) number. Your representative can also explain the available service replacement and repair programs.

The device should be packaged in its original transportation case and returned postage paid. A Smith & Nephew customer service representative will provide additional instructions for shipment.

Note: Product returned that is found to have been serviced by an unauthorized third-party repair facility and/or sterilized with a sterilization method other than one approved by Smith & Nephew will incur additional costs, regardless of warranty status.

Electrical Interference

This equipment is designed and tested to minimize interference with other electrical equipment. However, if interference occurs with other equipment, it may be corrected by one or more of the following measures:

- Reorient or relocate this equipment, the other equipment, or both.
- Increase the separation between the pieces of equipment.
- Connect the pieces of equipment into different outlets or circuits.
- Consult a biomedical engineer.

Environmental Protection

This equipment contains electronic printed circuit assemblies. At the end of the useful life of the equipment, it should be disposed of in accordance with any applicable national or institutional policy relating to obsolete electronic equipment.

Troubleshooting

Problem	Possible Cause	Suggested Action
Trauma Interface unit is without power	Mains power plug is not inserted (properly) or there is no mains power No power on the wall outlet One or both mains power fuses are blown	Insert mains power plug into reliable power supply Try another power outlet Replace mains fuses
Buttons or items are difficult to select on the touchscreen	Touchscreen is not calibrated	Access calibration software by selecting "Administration" from the "About" option under the "Menu" button. From the "Administration" screen, select the "Settings" button, then select the "Calibrate touch screen..." button
VGA video out not functioning	VGA port not activated on Trauma Interface Damaged VGA cable Video monitor not on correct input HDMI cable also connected	Connect VGA cable to both Trauma Interface and video monitor before powering on Trauma Interface Replace VGA cable Select proper input on video monitor Disconnect HDMI cable. Only one video output is possible at a time from the Trauma Interface
HDMI video out not functioning	HDMI port not activated on trauma interface Damaged HDMI cable Video monitor not on correct input VGA cable also connected	Connect HDMI cable to both Trauma Interface Replace HDMI cable Select proper input on video monitor Disconnect VGA cable. Only one video output is possible at a time from the Trauma Interface
TRIGEN™ SURESHOT™ Targeter not recognized (error message)	Error reading data from Targeter Connection issue with Targeter Damaged Targeter	Unplug Targeter, wait 10 seconds, and plug back in Power off Trauma Interface. Securely connect Targeter to unit then power on Trauma Interface Replace Targeter with new unit
Probe not recognized	Error reading data from probe Connection issue with probe Damaged probe	Unplug probe, wait 10 seconds, and plug back in Try other port on front of unit Power off Trauma Interface. Securely connect probe to unit then power on Trauma Interface Replace probe with new unit
Probe will not insert to the proper depth in the nail	Obstruction within the nail cannulation	Re-insert the ball tip guide rod into the nail cannulation to clear any obstruction

Problem	Possible Cause	Suggested Action
Nail not visible on the screen	Metal interference within the TRIGEN [®] SURESHOT [®] electromagnetic field TRIGEN SURESHOT Targeter and probe not within range of each other	Remove any metal objects from the targeting field Move the TRIGEN SURESHOT Targeter closer to the sensor end of the probe
Drill bit too short	Short drill bit being used and long drill sleeve option selected within software	Press "Menu", "Drill Sleeve" and select the short drill sleeve option to use the short drill bit
Drill bit too long	Long drill bit being used and short drill sleeve option selected within software	Press "Menu", "Drill Sleeve" and select the long drill sleeve option to use the long drill bit
Red and Green targeting circles representing the drill sleeve appear incorrect	Incorrect drill sleeve length selected Metal interference within the TRIGEN SURESHOT electromagnetic field Probe not inserted correctly within set stop Damaged probe	Verify the correct drill sleeve length is selected from the software menu Remove any metal objects from the targeting field Verify probe is oriented and seated correctly in the notches of the set stop Verify probe accuracy with Field Accuracy Gauge
Targeting missed the intended hole	Metal interference within the TRIGEN SURESHOT electromagnetic field Probe not inserted correctly within set stop Damaged probe	Remove any metal objects from the targeting field Verify probe is oriented and seated correctly in the notches of the set stop Verify probe accuracy with Field Accuracy Gauge
Drill sleeve cannot be removed from TRIGEN SURESHOT Targeter	Over-tightening of drill sleeve	Use the Slotted Hammer from the instrument tray as a wrench to unscrew the drill sleeve counter-clockwise from the Targeter
The 4.7 mm/4.0 mm step drill will not fit through the drill sleeve	Drill not compatible with the TRIGEN SURESHOT Distal Targeting System	Only use the long (7169-2811) and short (7169-2810) drills designed for use with the TRIGEN SURESHOT Distal Targeting System
Field accuracy check fails	Metal interference within the TRIGEN SURESHOT electromagnetic field Field Accuracy Gauge improperly installed on Targeter META Set Stop improperly installed on Field Accuracy Gauge Probe is incorrectly inserted within META Set Stop Probe is damaged Targeter is damaged	Remove any metal objects from the targeting field DO NOT perform test with Targeter still in the autoclave tray! Metal interference will occur Verify Field Accuracy Gauge is fully seated within Targeter port and knob is hand-tightened to Targeter Verify META Set Stop is correctly oriented and tightened securely to Field Accuracy Gauge Verify probe is oriented properly and inserted to the "REF" notches on probe body Replace probe with a new probe Replace Targeter with a new Targeter and return old one for service

Problem	Possible Cause	Suggested Action
Drill Depth Measurement will not trigger	Drill Depth Measurement Software is not enabled Operator may be trying to drill in a location that interferes with implant	Follow the instructions in the manual to enable the Drill Depth Measurement software Software will only allow Drill Depth Measurement activation to occur in areas that are safe to drill
Drill Depth Measurement triggers but cancels	Targeter is being rotated during action The angle to the nail hole is greater than 15 degrees	Do not rotate Targeter while performing Drill Depth Measurement feature Rotation is used to reset the Drill Depth Measurement Maintain alignment when targeting circles
Combo Drill Sleeve selection button is not available on the drill sleeve selection screen	Combo sleeve software feature is not enabled in the Administration menu	Follow the instructions under "Software configuration" in this User Manual to enable the Combo sleeve feature.
META-TAN™ nail selection button is not available on the implant selection screen	META-TAN nail software feature is not enabled in the Administration menu	Follow the instructions under "Software configuration" in this User Manual to enable the META-TAN Nail feature
On-screen Nail Rotation feature via the Targeter does not work	Nail Rotation software feature is not enabled in the Administration menu	Follow the instructions under "Software configuration" in this User Manual to enable the Nail Rotation feature
Unit does not display pictorial image of power switch when pressing the "power off" on-screen button. The screen simply turns black but the fan continues to run.	The "Shutdown V2" option has been disabled in the Administration menu	Follow the instructions under "Software configuration" in this User Manual to enable the "Shutdown V2" option Power off the unit via the rear power switch
All on-screen text is in a foreign language	A different language has been selected in the Administration menu	Follow the instructions under "Software configuration" in this User Manual to select the appropriate language from the available options

Technical Specifications

System power supply	Voltage Frequency Connected load	100–240 VAC 50–60 Hz ≤ 2.5 Amps
Fuses	Type Quantity	T3.15A H 250V 2
Classification	Protection class	Class I (with protective earth [=ground]) IPX0, continuous operation
	Type	BF (with hand-held FG, targeting probe)
	Humidity	No protection
	Explosion	No protection
Ambient conditions	Temperature	10–35°C
	Humidity	30–75% RH (non-condensing)
	Air pressure	700–1060 HPa
Storage and transport	Temperature	-20–50°C (in original packaging)
	Humidity	10–90% RH (non-condensing)
	Air pressure	700–1060 HPa
Interface	Dimensions (DxWxH)	31 cm x 26 cm x 13 cm ±10%
	Weight	3.5 kg ±10%
Transportation case	Dimensions (DxWxH)	276,35 x 355,6 x 269,75 mm±6,35 mm
Interface in Transportation case in corrugated transport shipper	Dimensions (DxWxH)	42 cm x 32 cm x 38 cm±10%
System + case	Weight	≈ 7 kg±10%
Complies with standards		USA, Canada, Europe IEC 60601-1 (IEC 60601-1:2005 Ed .3 + A1;C1:2014) IEC 60601-1-2 Ed. 4.0 (2014-02) CSA C22 .2#60601-1:2014 Ed .3

Electrical Safety

The system meets IEC standards (e.g. IEC 60601). All configurations comply with standard IEC 60601-1. Any person connecting equipment to the system is responsible for the configuration and must ensure that it complies with system standard 60601-1-1 or equivalent national standards. Please contact your customer service representative in the event of any queries.

WARNING: This equipment is not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide.

Cautions:

- The system may reboot after a power line surge greater than 1KV. This is considered a normal condition and the system will automatically reboot to the software application.
- Only use 110V/60Hz AC connection when powering the Trauma Interface within the United States. The Trauma Interface is not certified for 220V/60Hz operation within the United States.
- To ensure proper operation, no other electrical components should be located near the TRIGEN SURESHOT Distal Targeting unit.

Smith & Nephew reserves the right to make any technical changes.

Parts and Accessories List

Catalog No .	Description
71692851	TRIGEN SURESHOT Targeter 2nd Generation
71692852	Trauma Interface 2nd Generation
66800291	Power Cord, 250 Volt, 10 AMP, Continental Europe
71691166	SURESHOT Outer Screw Sleeve-Long
71691165	SURESHOT Inner Drill Sleeve-Long
71692804	SURESHOT Long Drill Sleeve
71692805	SURESHOT Short Drill Sleeve
71691154	SURESHOT Humeral 3.2 mm Drill Sleeve
71692814	SURESHOT META-NAIL Standard Drill Guide Probe
71692803	SURESHOT META-NAIL™ Semi-extended Drill Guide Probe
71692815	SURESHOT Percutaneous TAN/FAN Drill Guide Probe
71691152	SURESHOT Humeral Drill Guide Probe
71692806	SURESHOT META Set Stop
71692807	SURESHOT Percutaneous TAN™/FAN Set Stop
71691151	SURESHOT Humeral Set Stop
71692808	SURESHOT Field Accuracy Gauge
71692816	SURESHOT TAN Anteversion Locking Guide
71692809	SURESHOT Hexdriver
71691153	SURESHOT Humeral 3.5 mm Hexdriver
71691155	SURESHOT Humeral 3.2 mm AO Drill
71692811	SURESHOT Long AO Drill
71692810	SURESHOT Short AO Drill
71692832	TRIGEN SURESHOT Targeter Tray Insert 2nd Gen
71692830	TRIGEN SURESHOT Distal Targeting Instrument Tray
71692831	TRIGEN SURESHOT Distal Targeting Instrument Tray Lid
71691156	TRIGEN SURESHOT Humeral Instrument Caddy

CAUTION: Inspect all components regularly for wear. See Device Lifetime section of this document for additional information.

CAUTION: Use only Smith & Nephew disposables and accessories with the Smith & Nephew TRIGEN SURESHOT Distal Targeting System.

CAUTION: The TRIGEN SURESHOT probes are defined as single use items. The risk of reuse includes, but is not limited to, cross contamination between patients, probe malfunction and probe breakage. After use, the probe should be discarded per hospital policy/ procedures for biohazard material disposal.

Guidance and Manufacturer’s Declaration

Electromagnetic Emissions, Immunity and Separation Distances

Changes or modification to this system not expressly approved by the manufacturer may result in increased emissions or decreased immunity performance of the equipment or system and could cause EMC issues with this or other equipment. This system is designed and tested to comply with applicable regulations regarding EMC and shall be installed and put into service according to the EMC information stated as follows.

Cautions:

- Use of portable phones or other radio frequency (RF) emitting equipment near the system may cause unexpected or adverse operation.
- The presence of certain metal objects within the electromagnetic tracking volume system may cause unexpected or adverse operation.
- The equipment or system shall not be used adjacent to, or stacked with, other equipment. If adjacent or stacked use is necessary, the equipment or system shall be tested to verify normal operation in the configuration in which it is being used.
- The use of accessories, transducers and cables other than those specified may result in increased emissions or decreased immunity performance of the equipment or system.
- The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

Compliant Cables and Accessories

The table below lists cables, transducers, and other applicable accessories for which the manufacturer claims EMC compliance.

Note: Any supplied accessories that do not affect EMC compliance are not listed.

Part No .	Type	Description
7169-2851	Field Generator	TRIGEN SURESHOT Targeter 2nd Generation
7169-2814	Sensor	META-NAIL Standard Drill Guide Probe
7169-2815	Sensor	Percutaneous TAN/FAN Drill Guide Probe
7169-2803	Sensor	META-NAIL Semi-extended Drill Guide Probe
7169-1152	Sensor	Humeral Drill Guide Probe

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the TRIGEN SURESHOT Distal Targeting System, including cables specified by the manufacturer . Otherwise, degradation of the performance of this equipment could result .

Guidance and Manufacturer's Declaration – Electromagnetic Emissions

The TRIGEN SURESHOT Distal Targeting System is intended for use in the electromagnetic environment specified below. The customer or the user of the TRIGEN SURESHOT Distal Targeting System should assure that it is used in such an environment.

Emission Test	Compliance	Electromagnetic Environment - Guidance
RF Emissions CISPR 11	Group 1	The TRIGEN SURESHOT Distal Targeting System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The TRIGEN SURESHOT Distal Targeting System is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations / Flicker Emissions IEC 61000-3-3	Complies	

Guidance and Manufacturer's Declaration – Immunity

The TRIGEN SURESHOT Distal Targeting System is intended for use in the electromagnetic environment specified below. The customer or user of the TRIGEN SURESHOT Distal Targeting System should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
ESD IEC 61000-4-2	±8kV Contact ±2kV, ±4kV, ±8kV, ±15 kV air	±8kV Contact ±2kV, ±4kV, ±8kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the relative humidity should be at least 30%.
Electrical Fast Transient/Burst IEC 61000-4-4	±2kV 100kHz repetition frequency	±2kV 100kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5 kV, ±1 kV Line-to-line ±0.5 kV, ±1 kV, ±2 kV Line-to-Ground	±0.5 kV, ±1 kV Line-to-line ±0.5kV, ±1 kV, ±2kV Line-to-Ground	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the TRIGEN SURESHOT Distal Targeting System requires continued operation during power mains interruptions, it is recommended that the TRIGEN SURESHOT Distal Targeting System be powered by an uninterruptible power supply or battery.
Power frequency (50/60 Hz) Magnetic field IEC 61000-4-8	30A/M 50/60 Hz	30A/M 50/60 Hz	Power frequency magnetic fields should be that of a typical commercial or hospital environment.

Guidance and Manufacturer's Declaration – Immunity

The TRIGEN SURESHOT Distal Targeting System is intended for use in the electromagnetic environment specified below. The customer or user of the TRIGEN SURESHOT Distal Targeting System should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6Vrms (in ISM Bands) 150 kHz to 80 MHz	3Vrms 6Vrms (in ISM Bands) 150 kHz to 80 MHz	<p>Portable and mobile RF communications equipment should be used no closer to any part of the TRIGEN SURESHOT Distal Targeting System including cables than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ $d = 2.3 \sqrt{P}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters as determined by an electromagnetic site survey should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>
Radiated RF IEC 61000-4-3	80MHz to 2.7GHz	3 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz	

Guidance and Manufacturer's Declaration – Immunity						
Test Frequency (Mhz)	Band ^a (Mhz)	Service ^a	Modulation ^b	Maximum Power (W)	Distance (m)	Immunity Test Level (V/m)
385	380 - 390	TETRA 400	Pulse Modulation ^b 18 Hz	1.8	0.3	27
450	430 - 470	GMRS 460, FRS 460	FM ^c ± 5 kHz Deviation 1 kHz sine	2	0.3	28
710	704 - 787	LTE Band 13, 17	Pulse Modulation ^b 217 Hz	0.2	0.3	9
745						
780						
810	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse Modulation ^b 18 Hz	2	0.3	28
870						
930						
1 720	1 700 - 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse Modulation ^b 217 Hz	2	0.3	28
1 845						
1 970						
2 450	2 400 - 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse Modulation ^b 217 Hz	2	0.3	28
5 240	5 100 - 5 800	WLAN 802.11 a/n	Pulse Modulation ^b 217 Hz	0.2	0.3	9
5 500						
5 785						

Note: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the TRIGEN SURESHOT Distal Targeting System may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

- a. For some services, only the uplink frequencies are included.
- b. The carrier shall be modulated using a 50% duty cycle square wave signal.
- c. As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

Recommended separation distances between portable and mobile RF communications equipment and the TRIGEN SURESHOT Distal Targeting System

The TRIGEN SURESHOT Distal Targeting System is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the TRIGEN SURESHOT Distal Targeting System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the TRIGEN SURESHOT Distal Targeting System as recommended below, according to the maximum output power of the communication equipment.

Rated Maximum Output Power of Transmitter W	Separation Distance According to Frequency Transmitter m		
	150 kHz to 8 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Safe Disposal

All TRIGEN SURESHOT Distal Targeting System devices should be disposed or returned to Smith & Nephew. All devices associated with submitted complaints should be returned to Smith & Nephew for analysis. The TRIGEN SURESHOT Trauma Interface should be disposed of in accordance with any applicable national or institutional policy relating to obsolete electronic equipment. Reusable surgical devices that have been inspected and have reached the end of their lifetime should be sterilized, if possible, and disposed of according to institutional procedures. Devices that are unable to be sterilized before disposal should be disposed properly as a potential biohazard. Single-use devices should be disposed per institutional procedures for biohazard material disposal.

For Further Information

For further information, including warranty, please contact Customer Service at + 1 800 238 7538 for calls within the continental USA and +1 901 396 2121 for all international calls.

Smith & Nephew, Inc. 1450 Brooks
Road
Memphis, TN 38116 USA

+1 901 396 2121 International
+1 800 238 7538 U.S. Customer
Service

www.smith-nephew.com

ifu.smithnephew.com

©2017, 2020 Smith &
Nephew, Inc. All rights
reserved.
10/20 81103554 Rev.

ADDENDUM TO USER MANUAL
TRIGEN SURESHOT electromagnetic targeting surgical system for nails
distal locking during trauma operations with accessories

MD NAME, MODELS OR VARIATIONS OF MD

TRIGEN SURESHOT electromagnetic targeting surgical system for nails distal locking during trauma operations with accessories:

1. INTERFACE 2ND GENERATION in a protective soft case;
2. SURESHOT 2nd Gen TARGETER;
3. POWER CORD;
4. SURESHOT sleeves (if necessary):
 - TRIGEN SURESHOT Outer Screw Sleeve-Long;
 - TRIGEN SURESHOT Inner Drill Sleeve-Long;
 - SURESHOT Short Drill Sleeve;
 - SURESHOT Long Drill Sleeve;
 - TRIGEN SURESHOT Humeral 3.2 mm Drill Sleeve;
5. SURESHOT Guide Probes (starter kit):
 - SURESHOT META Probe - no more than 150 pcs.;
 - SURESHOT META SEMIEXTENDED Guide Probe - no more than 150 pcs.;
 - SURESHOT PERC TAN/FAN Probe - no more than 150 pcs.;
 - SURESHOT HUMERAL Probe - no more than 150 pcs..
6. SURESHOT AO Drills (starter kit):
 - SURESHOT Short AO Drill - no more than 150 pcs.;
 - SURESHOT Long AO Drill - no more than 150 pcs.;
 - TRIGEN SURESHOT Humeral 3.2 mm AO Drill - no more than 150 pcs..
7. SURESHOT Set Stops (if necessary):
 - SURESHOT Meta Set Stop;
 - SURESHOT Tan Set Stop;
 - SURESHOT Humeral Set Stop.
8. SURESHOT Field Accuracy Gauge (if necessary);
9. SURESHOT TAN Anteversion Locking Guide (if necessary);
10. SURESHOT Hexdriver (if necessary):
 - SURESHOT Hexdriver 4.7mm;
 - TRIGEN SURESHOT Humeral 3.5 mm Hexdriver.

Accessories:

1. TRIGEN SURESHOT Targeter Tray Insert 2nd Gen;
2. Trigen SURESHOT Dist Target Instr Tray;
3. Trigen SURESHOT Dist Target Inst Tr Lid;
4. TRIGEN SURESHOT Humeral Instrument Caddy.

ADVERSE EFFECTS

There are no potential adverse effects.

THE POSSIBILITY AND WAYS OF INTEGRATING WITH OTHER MEDICAL DEVICES

The SURESHOT system is intended for use together with Smith & Nephew implants, intramedullary nails and screws.

DESCRIPTION OF CLEANING, DISINFECTION AND STERILIZATION PROCESS

Non-sterile reusable components and accessories (Sleeves, Field Accuracy Gauge, targeter, Anteversion Locking Guide, Set Stops, Hexdrivers)

Non-sterile reusable system components that need to be sterilized fall into the following product groups as defined in document «Instructions for Care, Maintenance, Cleaning and Sterilization of Smith & Nephew Orthopaedics Devices»:

Products	Cleaning Product Groups
TRIGEN SURESHOT Outer Screw Sleeve-Long TRIGEN SURESHOT Inner Drill Sleeve-Long SURESHOT Short Drill Sleeve SURESHOT Long Drill Sleeve TRIGEN SURESHOT Humeral 3.2 mm Drill Sleeve SURESHOT Meta Set Stop	Devices with challenging design features

SURESHOT Tan Set Stop SURESHOT Humeral Set Stop SURESHOT Field Accuracy Gauge SURESHOT TAN Anteversion Locking Guide SURESHOT Hexdriver 4.7mm TRIGEN SURESHOT Humeral 3.5 mm Hexdriver	
SURESHOT 2nd Gen TARGETER	Powered devices

Accessories that are used to sterilize components as ancillary items belong to the following product groups, as defined in document «Instructions for Care, Maintenance, Cleaning and Sterilization of Smith & Nephew Orthopaedics Devices»:

Products	Cleaning Product Groups
TRIGEN SURESHOT Targeter Tray Insert 2nd Gen Trigen SURESHOT Dist Target Instr Tray Trigen SURESHOT Dist Target Inst Tr Lid TRIGEN SURESHOT Humeral Instrument Caddy	Containment devices

Refer to the document «Instructions for Care, Maintenance, Cleaning and Sterilization of Smith & Nephew Orthopaedics Devices» for procedures for cleaning, disinfecting and sterilizing

Sterile disposable components (Probes, Drills)

	STERILIZATION METHOD FOR SINGLE USE COMPONENTS PROVIDED STERILE
SURESHOT META probe SURESHOT SEMI-EXTENDED guide probe SURESHOT perc TAN/FAN probe SURESHOT HUMERAL probe SURESHOT Humeral Drill Guide Probe	Ethylene oxide
SURESHOT Short AO Drill SURESHOT Long AO Drill SURESHOT Humeral 3 2 mm AO Drill	Gamma ray

HOW SUPPLIED

All components are supplied in 1 piece per package.

The TRIGEN SURESHOT Interface 2nd generation in a soft case and intended powerd cord are supplied separately.

The TRIGEN SURESHOT Targeter Tray Insert 2nd Gen comes with 3 screws, a screwdriver for installation and a sticker.

A Russian-language user manual in paper version are supplied if necessary. A Russian-language user manual is available on the website by the link indicated in the Russian-language labelling.

The TRIGEN SURESHOT Distal Targeting System devices are packaged separately, with the exception of some reusable surgical devices that may be provided as a set within an instrument tray

LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN

The subject Smith & Nephew products do not incorporate non-viable materials of animal origin, nor are they present in the manufacturing process or present in the final device; therefore this section does not apply.

INFORMATION ABOUT MEDICINES AND (OR) PHARMACEUTICAL SUBSTANCES INCLUDED IN MD

The subject Smith & Nephew products do not incorporate a medicinal substance, nor is it present in the manufacturing process or present in the final device; therefore this section does not apply.

OPERATION CONDITIONS FOR WORKING PARTS THAT COME IN CONTACT WITH THE INTERNAL ENVIRONMENT OF THE BODY

Working parts that come in contact with the internal environment of the body can be operated in the temperature range 32-42 °C.

Standard	Description
BS EN 62366:2008	Medical Devices. Application of usability engineering to medical devices
IEC 60601-1-6:2010	Medical Electrical Equipment –Part 1-6: General Requirements for basic safety and essential performance- Collateral standard: Usability
BS EN ISO 14971: 2012	Medical devices. Application of risk management to medical devices
BS EN ISO 16061:2015	Instrumentation for use in association with non- active surgical implants - General requirements
BS EN ISO 13485:2016	Medical devices. Quality management systems. Requirements for regulatory purposes
IEC 60601-1:2005 +A1: 2012	Medical electrical equipment-Part 1: General Requirements for basic safety and essential performance
BS EN ISO 11607-2:2006+A1:2014	Packaging for terminally sterilized medical devices. Validation requirements for forming, sealing and assembly processes
BS EN ISO 11607-2:2006+A1:2014	Packaging for terminally sterilized medical devices. Validation requirements for forming, sealing and assembly processes
ICH GCP E6, E8, E9, E10	Good Clinical Practice
BS EN ISO 14155:2011	Clinical investigation of medical devices for human subjects - Good clinical practice
NB-MED 2.12-1 rev. 8	Post-Marketing Surveillance (PMS) post market/production
MEDDEV 2.12/2 rev. 2	Post Market Clinical Follow-up Studies
MEDDEV 2.7.1/Rev. 4	Clinical Evaluation
ISO 11137-1:2006	<i>Sterilization of health care products - Requirements for validation and routine control radiation sterilization</i>
ISO 11135-1:2007	Sterilization of health care products - Ethylene oxide - Part 1: Requirements for development, validation, and routine control of a sterilization process for medical devices.
ISO 10993-1: 2009	Biological Evaluation of Medical Devices - Part 1: Evaluation and testing

WARRANTY AND RECLAMATION

Smith & Nephew LLC (S&N) warrants all of its products to the original purchaser (Customer) against defect in workmanship and materials for one hundred eighty (180) days from date of purchase from Smith & Nephew; except for NAVIO® products, one (1) year from date of installation (Warranty Period). For any products found to not be in conformance with this warranty during the Warranty Period, this warranty provides and is restricted to, as elected by S&N, either (i) repair or replacement of such products without charge and within a reasonable period of time or (ii) a refund or credit in the amount of the purchase price of such products. If S&N repairs or replaces product under this warranty and requests Customer to return such product, Customer must ship such product to S&N freight prepaid by Customer. Customer shall be invoiced for any replacement product if Customer does not return the requested replaced product within thirty (30) days after S&N's shipment of the replacement product. This warranty does not cover and is voided by any of the following: (i) product packaged or labeled by someone other than S&N or its authorized agents; (ii) product not used in compliance with the specifications, instructions or claims for use of the product; (iii) equipment product used in conjunction with disposables or accessories not specified for use with such equipment; (iv) equipment

product used in conjunction with reprocessed disposables or accessories; (v) modification of product;

(vi) product past its expiration date; (vii) normal wear and tear; (viii) damage due to misuse, reprocessing, alteration, unauthorized repair or negligent handling or damage due to lack of care the owner, user or handler of the product including but not limited to storage, handling or cleaning; or (ix) any other damage inflicted to products by the owner, user or handler. A 90-day warranty on repairs applies to the defective component repaired. This warranty applies only to the original purchaser from S&N (or its authorized distributor) and is not transferable. EXCEPT TO THE EXTENT PROHIBITED OR OTHERWISE REQUIRED BY APPLICABLE LAW, THIS WARRANTY IS THE SOLE WARRANTY OF S&N, AND ALL OTHER WARRANTIES OF ANY KIND OR DESCRIPTION WHATSOEVER, INCLUDING WARRANTIES OF MERCHANTABILITY, SATISFACTORY QUALITY AND FITNESS FOR A PARTICULAR PURPOSE, EXPRESSED OR IMPLIED, ARE EXCLUDED.

AUTHORIZED REPRESENTATIVE ON THE TERRITORY OF THE RUSSIAN FEDERATION: Smith & Nephew, LLC, Russia, 105120, Moscow, 2nd Syromyatnichesky lane., 1, floor 9, prem. 1, room 1 Tel: +7 (495) 755-55-03

NAME AND LEGAL ADDRESS OF MANUFACTURER AND AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION

DEVELOPER:

Smith & Nephew, Inc.

1450 Brooks Road Memphis, Tennessee (TN) 38116, USA (CWA)

MANUFACTURER:

Smith & Nephew, Inc.

1450 Brooks Road Memphis, Tennessee (TN) 38116, USA (CWA)

AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION

Smith & Nephew, LLC, Russia, 105120, Moscow, 2nd Syromyatnichesky lane., 1, floor 9, prem. 1, room 1

Tel: +7 (495) 755-55-03

Instructions for Care,
Maintenance, Cleaning and
Sterilization of Smith & Nephew
Orthopaedics Devices


Table of Contents

Purpose	4	Containment devices	10
Scope	4	Automatic washing cycle steps and parameters	10
Definitions	4	Disinfection	10
Recommended Cleaning Instructions	5	Verifying cleaning	10
Important Information/Recommendations for Use	5	Inspection and function testing	10
Warnings	6	Maintenance and Care	11
Limitations on Processing	6	Storage Between Cleaning and Sterilization	11
Limitations of cleaning instructions	6	Recommended Sterilization Instructions	11
Limitations of processing	6	Warnings	11
Cleaning Product Groups	6	Limitations of Processing	11
Overview of product groups for reusable device cleaning	6	Implants	11
Definition of product groups	6	Containment Devices	11
Initial Treatment at the Point of Use	7	Preparation for Sterilization	12
Transportation	7	Reusable Medical Devices	12
Preparation Before Cleaning	7	Single Use Medical Devices	12
Water Quality for Cleaning	7	Sterilization Wrap/Reusable Rigid Containers	12
Cleaning: Manual	7	Water Quality	12
Devices Without challenging design features	7	Recommended Sterilization Parameters	13
Devices With challenging design features	7	Dynamic Air Removal (Prevacuum) Steam	13
Flexible reamer devices	8	Gravity Displacement Steam	13
Powered Devices	8	For Non-US Customers - UK Steam Cycle Prevacuum cycle	13
Containment Devices	8	World Health Organization (WHO) Steam Cycle	13
Non-critical Devices	8	Storage After Cleaning and Sterilization	13
Verifying Cleaning	9	Cleaning and Sterilization Validation Information	14
Drying	9	References	14
Cleaning: Automated	9	Contact Information	14
Devices Without challenging design features	9	Glossary of Symbols	16
Devices With challenging design features	9		
Flexible reamer devices	9		
Powered devices	10		

Purpose/Scope/Definitions

Purpose

This document was prepared to provide instructions for the care, maintenance, cleaning, and sterilization of the medical devices produced by the Orthopaedics Business of Smith & Nephew.

These methods were developed using standard equipment and practices common to health care facilities. Validation testing to support these instructions was based on recognized guidelines and standards for reusable devices and containment devices from the following organizations:

- American National Standards Institute (ANSI)
- Association for the Advancement of Medical Instrumentation (AAMI)
- Association of Operating Room Nurses (AORN)
- German Instrument Working Group (AKI) Arbeitskreis Instrumenten-Aufbereitung
- Health Canada
- International Standards Organization (ISO)
- International Association of Healthcare Central Service Material Management (IAHCSMM)
- World Health Organization (WHO)
- UK Department of Health
- Australian/New Zealand Standard
- Centers for Disease Control (CDC)
- U.S. Food and Drug Administration (FDA)

Scope

These instructions apply to all reusable medical devices that are sold by Smith & Nephew Orthopaedics that are initially supplied nonsterile or sterile and require the user to process them after initial/subsequent use. These instructions also apply to single use medical devices that are supplied sterile and nonsterile and require processing. Both sterile and nonsterile single use medical devices i.e. plates, nails, screws, pins and wires are commonly placed in containment devices for use and therefore require processing prior to use.

These instructions apply to Smith & Nephew Orthopaedic devices that fall under Spaulding's classification scheme for medical devices as "Critical Devices" and "Noncritical Devices". "Critical Devices" are introduced directly into the human body either into or in contact with the bloodstream or into other normally sterile areas of the body. "Critical Devices" present a high degree of risk of transmission of infection if contaminated and therefore must be sterile at the time of use. "Noncritical Devices" contact only the intact skin of the patient.

Definitions

Biological Indicator (BI): Test system that contains viable microorganisms that provide a defined resistance to a specified sterilization process.

Caddy: A small flip top containment device that contains multiple small implants (i.e. screws, plates, etc.).

Chemical Indicator (CI): Test system that reveals change in one or more pre-defined process variables based on a chemical or physical change resulting from exposure to a process. CIs assist in the detection of potential sterilization failures that could result from incorrect packaging, incorrect loading of the sterilizer, or sterilizer malfunctions. The "pass" response of a CI does not prove that the item monitored by the indicator is sterile.

Cleaning: The removal of contaminants to the extent necessary for further processing or for intended use.

Containment Device: A reusable rigid sterilization container, instrument case/cassette, caddy, or organizing tray and any reusable accessories for use in healthcare facilities for the purpose of containing reusable medical devices for sterilization.

Critical Water: Water that has been extensively treated (usually by a multi-step treatment process that could include a carbon bed, softening, Deionized (DI), and Reverse Osmosis (RO) or distillation) to ensure that the microorganisms and the inorganic and organic materials are removed from the water. A final submicron filtration could also be part of the treatment process. This water is mainly used for the final rinse or steam generation.

Disinfection: A process to reduce the number of viable microorganisms to a level previously specified as being appropriate for a defined purpose.

Immediate Use: The shortest possible time between a sterilized item's removal from the sterilizer and its aseptic transfer to the sterile field. A sterilized item intended for immediate use is not stored for future use or held from one case to another.

Manual Cleaning: The removal of contaminants from an item to the extent necessary for further processing or for intended use without the use of an automated process. Ultrasonic cleaning is considered a manual cleaning step for the purposes of this document.

Manufacturer: The natural or legal person with responsibility for design and/or manufacture of a medical device with the intention of making the medical device available for use, under his name; whether or not such a medical device is designed and/or manufactured by that person himself or on his behalf of another person(s).

Medical Device: Instrument, apparatus, implant, reagent for in-vitro use, software, material or other similar or related article intended by the manufacturer to be used alone or in combination, for human beings for one or more of the specific medical purpose(s) of:

- diagnosis, prevention, monitoring, treatment or alleviation of disease; diagnosis, monitoring, treatment, alleviation of or compensation for an injury; investigation, replacement, modification or support of the anatomy or of a physiological process; supporting or sustaining of life; control of conception; disinfection of medical devices; providing information by means of in vitro examination of specimens derived from the human body; and does not achieve its primary intended action by pharmacological, immunological or metabolic means, in or on the human body, but which may be assisted in its function by such means.

Processing: Preparation of medical devices or activity to prepare a new or used healthcare products for its intended use. In this document processing includes cleaning, disinfection and sterilization.

Reusable Medical Device: A medical device designated or intended by the manufacturer as suitable for processing and reuse.

Single Use Medical Device (SUD): A medical device that is designated or intended by the manufacturer for one-time use only. For clarity, a single use device that has come into contact with blood, tissue, or bodily fluids is not intended to be further processed and used again.

Sterile: Free from viable microorganisms.

Sterility Assurance Level (SAL): The probability of a single viable microorganism occurring on an item after sterilization, expressed as the negative exponent to the base 10.

Sterilization: The process used to render product free from viable microorganisms.

Utility Water: Water as it comes from the tap that might require further treatment to achieve the specifications. This water is mainly used for flushing, washing, and rinsing.

Validation: Confirmation, through the provision of objective evidence, that the requirements for a specific intended use or application have been fulfilled.

Washer-Disinfector: Equipment that is designed to clean and disinfect product.

Recommended Cleaning Instructions

Cleaning is the single most important step in preparing a device for reuse. Effective cleaning must be carried out to achieve proper disinfection/sterilization. Cleaning begins at the point of use to prevent soil and contaminants from drying on the medical devices after use. Thorough cleaning and rinsing are vital to processing reusable medical devices. Also, thorough rinsing is important for the removal of any residual cleaning agents from the medical devices. The purpose of cleaning and rinsing is to remove all adherent visible debris and to reduce the number of particulates, microorganisms, and pyrogens. The recommended cleaning instructions in this document include both manual and automatic washing/disinfection procedures. While manual cleaning is the most universal method of cleaning, automatic washing is preferred. The cleaning processes presented in this brochure have been validated. Other methods of cleaning may be suitable but must be validated by the user of the device.

Important Information/Recommendations for Use

Automatic washer/disinfector

Washer-disinfectors are not only used to clean devices, but also to provide intermediate to high level disinfection with a hot water rinse. Cleaning is dependent upon thorough coverage of the devices and the force of the water spray. Therefore, all sections of the device must be accessible for ease of cleaning and penetration of cleaning agents. The automatic washer/disinfector equipment should be operated following the manufacturer's instructions for use.

Detergents

Low foaming detergents with a pH range between 6.0 and 8.0 are recommended. Detergents with a pH outside this range can have an adverse effect or be damaging to some medical devices and containment devices. Enzymatic detergents aid in the removal of organic soil such as blood.

Detergents should be used at the concentration and temperature recommended by the detergent manufacturer.

Some alkaline detergents have been formulated to be safe for processing medical devices. Refer to the manufacturer's limitations and warnings for information concerning specific materials that are adversely affected by the detergents.

Loaner Instrument Set

It is the hospital's responsibility to clean, disinfect, package and sterilize all loaner devices upon receipt. In addition, it is also the responsibility of the hospital to clean, disinfect, package and sterilize the instrument sets prior to their return to Smith & Nephew.

Recommended Cleaning Instructions *continued*

Manual cleaning tools

Hospital tools necessary for manual cleaning include: surgical scrub brushes, chenille pipe cleaners, soft low linting cloths, cotton tip applicators, and several various sizes and lengths of brushes.

Do not use abrasive cleaning tools (i.e. scouring pads or metal brushes).

Cleaning tools must be cleaned and inspected after use. Cloths should be clean and lint free and changed frequently. Brushes should be clean. Discard worn brushes and disposable cleaning tools.

Safety precautions

Personal Protective Equipment (PPE) should be worn when handling or working with contaminated or potentially contaminated devices. PPE includes: gown, mask, goggles or face shield, and shoe covers.

Universal precautions are standards of infection control practices designed to reduce the risk of transmission of bloodborne infections. Universal precautions should be observed by all hospital personnel that work with contaminated or potentially contaminated devices. Exercise caution when handling devices with sharp points or cutting edges.

Water

The quality of water should be carefully considered for use in cleaning reusable devices. The water quality can affect the life of the device. Water hardness is a concern because deposits left on medical devices may result in ineffective cleaning and sterilization. Final rinsing should be carried out using critical water.

Warnings

- It is the responsibility of the user to ensure that the cleaning process is performed following these procedures to achieve the desired result.
- For ultrasonic cleaning, the enzymatic detergent solution should be changed before it becomes heavily soiled so that effective cleaning is not inhibited.
- Most electronic devices cannot be submerged during cleaning and could sustain permanent damage as a result of submersion. Check the "Instructions for Use" for information on submersion of electronics.
- Do not sonicate powered devices.
- All cleaning should be performed in a manner designed to minimize exposure to bloodborne pathogens. Manual cleaning should be done while the instrument is immersed.

Limitations on Processing

Limitations of cleaning instructions

These recommended procedures are intended as a general guide for cleaning of medical devices. Some devices are labeled with more specific instructions. Refer to the Instruction for Use for that device.

Limitations of processing

Repeated processing of reusable medical devices has minimal effect on the devices.

The useful life of these devices depends on many factors including the method and duration of each use, and the handling between uses. Careful inspection and functional test of the device before use is the best method of determining the end of serviceable life for the medical device. Do not use devices that show evidence of damage and wear. Evidence of damage and wear on a device may include but is not limited to corrosion (i.e. rust, pitting), discoloration, excessive scratches, flaking, wear and cracks. Improperly functioning devices, devices with unrecognizable markings, missing or removed (buffed off) part numbers, damaged and excessively worn devices should not be used.

Cleaning Product Groups

Overview of product groups for reusable device cleaning

Cleaning of reusable devices is dependent upon product design features. The cleaning methods for Smith & Nephew Orthopaedics reusable devices are based on product groups that have design features that present a similar challenge to cleaning.

Cleaning Product Groups:


- Devices **without** challenging design features
- Devices **with** challenging design features
- Flexible reamers
- Powered devices
- Containment devices
- Noncritical devices

Definition of product groups

Devices without challenging design features: Includes all devices that do not have design features that present a challenge to cleaning by the Smith & Nephew Orthopaedics recommended cleaning procedure. These devices do not have difficult to access locations for cleaning such as lumens, interfaces, hinged/mating surfaces, crevices, holes, serrations, etc. These devices do not have retractable or moving parts.

Recommended Cleaning Instructions *continued*

Examples: Bone spikes, osteotomes, mallets.

Devices with challenging design features: Includes all devices that have design features that present a challenge to cleaning by Smith & Nephew Orthopaedics recommended cleaning procedure such as lumens, interfaces, hinged/mating surfaces, crevices, holes, and serrations etc. These devices may have retractable and moving features.

Examples: Reamers, T-handles, cable passers, cutting blocks, and hinged clamps.

Flexible reamers: Includes all devices that have a laser cut flexible reamer design.

Examples: Flexible screw drills and flexible shafts.

Powered devices: Those that use pneumatic power or have power cords and require the use of electricity to operate.

Examples: ACCURIS® handpiece and pulse lavage handpiece.

Containment devices: Instrument cases, trays, caddies and lids

Non-critical devices: Those that do not contact the patient or may contact intact skin.

Examples: Power supply boxes and power cable surfaces.

Initial Treatment at the Point of Use

Cleaning begins at the point of use to prevent soil and contaminants from drying on the devices. Keep the devices moistened after use and before cleaning.

Transportation

Used devices must be transported to the processing location (i.e. Sterile Processing) in closed or covered containers to prevent contamination risks.

Preparation Before Cleaning

Disassemble each device into its component parts. Containment devices must be cleaned separately from the medical devices. If you have questions concerning the disassembly of any Smith & Nephew Orthopaedics device, contact your Smith & Nephew sales representative.

Water Quality for Cleaning

Use utility water for flushing, washing and rinsing during cleaning except for the final rinse. For the cleaning final rinse use critical water.

Cleaning: Manual

Devices Without challenging design features

1. Rinse in cold water <45° C (113° F) to remove visible debris and to prevent coagulation of blood.
2. Immerse and soak for a minimum of one (1) minute in enzymatic detergent.
3. With the device completely immersed in the cleaning solution, use a surgical scrub brush to remove visible debris from the surfaces of the device.
4. Rinse for a minimum of one minute using a sufficient volume of critical water that will completely immerse the device. Water should be changed a minimum of two times to ensure thorough rinsing.
5. Check for visible debris. Repeat cleaning if debris is visible.

Devices With challenging design features

1. Rinse in cold water <45° C (113° F) to remove visible debris and to prevent coagulation of blood.
2. Immerse and soak for a minimum of five (5) minutes in enzymatic detergent.
3. With the device completely immersed in the cleaning solution, remove additional soil from the surfaces and from the challenging design features (i.e. holes, lumens, hinged/mating surfaces, interfaces, crevices, serrations) using common hospital cleaning tools.
 - a. Move and/or retract all moveable parts and remove visible debris using a brush.

Recommended Cleaning Instructions *continued*

- b. Scrub lumens or holes with a brush of an appropriate size to ensure that the full width and depth is accessed. Use a twisting action with the brush. Thoroughly flush lumens with enzymatic detergent using a syringe.
 - c. Open hinged devices and scrub hinged area with a brush.
 - d. Scrub crevices with a brush.
 - e. Scrub the surfaces of the device with a toothbrush styled brush.
4. Sonicate the device in enzymatic detergent in its fully opened position for a minimum of 15 minutes in an ultrasonic cleaner.
 5. Rinse for a minimum of one minute using a sufficient volume of critical water that will completely immerse the device. Water should be changed a minimum of two times to ensure thorough rinsing. If the components of the device are movable or can be retracted, it is necessary to do so during the rinsing process. Thoroughly flush lumens with critical water using a syringe.
 6. Check devices for visible debris (see "Verifying Cleaning"). Repeat cleaning if debris is visible.

Flexible reamer devices

1. Gently bend the flexible reamer in a U-shape and rinse in cold water <45° C (113° F) to remove visible debris and to prevent coagulation of blood.
2. Immerse and soak for a minimum of ten (10) minutes in enzymatic detergent.
3. With the device completely immersed in the cleaning solution, scrub lumen with a brush of an appropriate size to ensure that the full width and depth is accessed. Use a twisting action with the brush. (The bristles should be stiff enough to remove bone and tissue.) Thoroughly flush lumens with enzymatic detergent using a syringe.
4. With the device completely immersed in the cleaning solution, scrub the surface with a surgical scrub brush to remove all visible debris from the surface and crevices.
5. With the device completely immersed in the cleaning solution, gently bend the reamer in a U-shape and scrub the surface with a scrub brush.
Note: Gently bend at several locations along the length to access all crevices.
6. Gently bend the reamer in a U-shape and rinse thoroughly with warm water making sure to irrigate the lumens and crevices. Thoroughly flush lumens with water using a syringe.
Note: Use a clean brush during the rinse cycle and move it back and forth through the lumen several times during rinsing.
7. Sonicate in enzymatic detergent for a minimum of 15 minutes in an ultrasonic cleaner.

8. Rinse for a minimum of one minute using a sufficient volume of critical water to completely immerse the instrument. Water should be changed a minimum of two times to ensure thorough rinsing. Thoroughly flush lumens with critical water using a syringe.
9. Check the devices for visible debris (see "Verifying Cleaning"). If debris is visible, repeat cleaning steps 2 – 8.

Powered Devices

Leave the device that directly connects to the motorized device (i.e. the hose that connects to a powered handpiece) in place during cleaning to prevent an excessive amount of water from entering the motor.

Do not sonicate powered devices.

Clean the Motorized Device as Follows:

1. Rinse in cold water <45° C (113° F) to remove visible debris and to prevent coagulation of blood.
2. With the connecting device still in place, immerse the motorized device, with the connecting instrument still in place, in enzymatic detergent and soak for three minutes.
3. With the device completely immersed in the cleaning solution, scrub the surfaces of the device in the enzymatic detergent using cleaning brushes.
4. Rinse for a minimum of one minute using a sufficient volume of critical water to completely immerse the instrument. Water should be changed a minimum of two times to ensure thorough rinsing.
5. Remove the motorized device from the water before disassembly of the connecting device.
Clean the connecting device separately following the cleaning instructions for the correct Cleaning Product Group.

Containment Devices

1. Remove medical devices from the containment device prior to cleaning. Inspect the containment device for visible debris.
2. If dried visible debris is observed, follow the manual instructions for "Devices **with** challenging design features."
3. If dried visible debris is not observed, clean following the manual instructions for "Devices **without** challenging design features."

Non-critical Devices

Sterile hospital grade covers can be used whenever possible to prevent contamination of non-critical devices. Porous items like foam that are contaminated must be discarded. Clean surfaces by wiping with a cloth and using detergent and water.

Verifying Cleaning

1. After cleaning, visually inspect devices under normal lighting for any remaining visible debris.
2. For difficult to view design features, apply 3% hydrogen peroxide. Bubbling indicates the presence of blood. Rinse the devices thoroughly with critical water following hydrogen peroxide testing.
3. If not visibly clean, repeat cleaning and reinspect the device.

Drying

Use clean filtered compressed air or a clean lint free towel to remove moisture from the devices.

Cleaning: Automated

Devices *Without* challenging design features

1. Manual precleaning:
IS NOT REQUIRED if the device does not have dried-on visible debris. Place the device directly into the automatic washer for cleaning.
IS REQUIRED if the device does have dried on visible debris. Follow the manual cleaning steps below prior to placing the device in the automatic washer for cleaning.
 - a. Rinse in cold water <45° C (113° F) to remove visible debris and to prevent coagulation of blood.
 - b. Immerse and soak for a minimum of one (1) minute in enzymatic detergent.
 - c. With the device completely immersed in the cleaning solution, use a surgical scrub brush to remove visible debris.
 - d. Rinse for a minimum of one minute using a sufficient volume of water to completely immerse the device. Water should be changed a minimum of two times to ensure thorough rinsing.
2. Place the device in the automatic washer and run the recommended automatic washer steps (see section "Automatic Washing Cycle Steps and Parameters.")

Devices *With* challenging design features

1. Manual precleaning is required for all devices in this product group. Follow the manual cleaning steps below prior to placing the device in the automatic washer for cleaning.
 - a. Rinse in cold water <45° C (113° F) to remove visible debris and to prevent coagulation of blood.
 - b. Immerse and soak for a minimum of five (5) minutes in enzymatic detergent.

- c. With the device completely immersed in the cleaning solution, remove additional visible debris from the device surfaces and from the challenging design features (i.e. holes, lumens, hinged/mating surfaces, interfaces, crevices, serrations) using common hospital cleaning tools.
 - Move and/or retract all moveable parts and use a brush to remove visible debris.
 - Scrub lumens and holes with a brush of an appropriate size to ensure that the full width and depth is accessed. Use a twisting action with the brush. Thoroughly flush lumens with enzymatic detergent using a syringe.
 - Open hinged devices and scrub hinged area with a brush.
 - Scrub crevices with a brush.
 - d. Sonicate the device in its fully opened position for a minimum of 15 minutes in an ultrasonic cleaner containing enzymatic detergent.
 - e. Rinse for a minimum of one minute using a sufficient volume of water to completely immerse the device. If the components of the instrument are movable or can be retracted, do so during the rinsing process. Change the water a minimum of two times to ensure thorough rinsing. Thoroughly flush lumens with water using a syringe.
2. Load the devices into the washer such that all design features of the device are accessible to cleaning and such that design features that might retain liquid can drain (for example, hinges should be open and lumens and holes positioned to drain).
 3. Run the recommended automatic washer steps (see section "Automatic Washing Cycle Steps and Parameters.")
 4. Check devices for visible debris (see "Verifying cleaning"). Repeat cleaning if debris is visible and reinspect.

Flexible reamer devices

1. Manual precleaning is required for all devices in this product group. Follow the manual cleaning steps below prior to placing the device in the automatic washer for cleaning.
 - a. Gently bend the flexible reamer in a U-shape and rinse in cold water <45° C (113° F) to remove visible debris and to prevent coagulation of blood.
 - b. Immerse and soak for a minimum of ten (10) minutes in enzymatic detergent.
 - c. With the device completely immersed in the cleaning solution, scrub the lumen with a brush of an appropriate size to ensure that the full width and depth is accessed. Use a twisting action with the brush. (The bristles should be stiff enough to remove bone and tissue.)
 - d. With the device completely immersed in the cleaning solution, scrub the surface with a surgical scrub brush to remove all visible debris from the surface and crevices. Thoroughly flush lumens with enzymatic detergent using a syringe.

Recommended Cleaning Instructions *continued*

- e. With the device completely immersed, gently bend the reamer in a U-shape and scrub the surface with a scrub brush. **Note:** Bend at several locations along the length to access all crevices.
 - f. Gently bend the reamer in a U-shape and rinse for a minimum of one minute using a sufficient volume of water to completely immerse the device. Change the water a minimum of two times to ensure thorough rinsing. Make sure to irrigate the lumen and crevices. Thoroughly flush lumens with enzymatic detergent using a syringe.
 - g. Sonicate for a minimum of 15 minutes in an ultrasonic cleaner containing enzymatic detergent.
2. Place the device in the automatic washer and run the recommended automatic washer steps (see the section "Automatic Washing Cycle Steps and Parameters").

Powered devices

Leave the device that directly connects to the motorized device (i.e. the hose that connects to a powered handpiece) in place during cleaning to prevent an excessive amount of water from entering the motor.

Do not sonicate powered device.

Clean the motorized device as follows:

1. Rinse in cold water <45° C (113° F) to remove visible debris and to prevent coagulation of blood.
2. With the device completely immersed in the cleaning solution, use a surgical scrub brush to remove visible soil from the surfaces of the device and connecting device.
3. Rinse for a minimum of one minute using a sufficient volume of water to completely immerse the device. Change the water a minimum of two times to ensure thorough rinsing.
4. Place the motorized device with the connecting device attached to it in the automatic washer and run the recommended automatic washing steps (see the section "Automatic Washing Cycle Steps and Parameters").

Clean the connecting device separately following the cleaning instructions for the correct Cleaning Product Group.

Containment devices

1. Remove all medical devices from the containment device and inspect it for visible debris.
2. If dried visible debris is observed follow the manual cleaning instructions for "Devices with challenging design features".
 - a. Then place the containment device in the automatic washer and run the recommended automatic washer steps (see the section "Automatic Washing Cycle Steps and Parameters").
3. If dried debris is not observed, place the containment device in the automatic washer and run the automatic washer steps given below (see section "Automatic Washing Cycle Steps and Parameters").

Automatic washing cycle steps and parameters

Load the devices into the automatic washer-disinfector.

- Connect any device lumen(s) to the rinsing ports of the washer-disinfector. If no direct connection is possible, place the lumen(s) on injector jets or in injector sleeves of the injector basket.
- Place blind holes at a downwards incline for drainage.
- Place all instruments in the open position whenever possible (i.e. hinged/mating features, articulating surfaces).

Minimum cycle parameters:

- Five (5) minute cold prewash <45° C (113° F)
- Five (5) minute neutral enzymatic detergent wash
- Five (5) minute neutral detergent wash
- One (1) minute rinse with critical water.

Disinfection

All Smith & Nephew orthopaedics reusable devices can be thermally disinfected using the parameters below.

- minimum 91° C for 1 minute
- 25 minute minimum dry time

Use clean filtered compressed air or a clean lint free towel to remove remaining moisture from the devices.

Verifying cleaning

1. After cleaning, visually inspect devices under normal lighting for the removal of visible debris.
2. For difficult to view design features, apply 3% hydrogen peroxide. Bubbling indicates the presence of blood. Rinse the devices thoroughly with critical water following hydrogen peroxide testing.
3. Repeat cleaning if not visibly clean and reinspect the device.

Inspection and function testing

All reusable devices: Visually inspect for damage or wear, including components in their disassembled state prior to re-assembly. Ensure components are re-assembled securely.

Hinged devices: Check for smooth movement of the hinge without excessive "play."

Locking mechanisms: Check for action.

Cutting features: Check edges for distortion and/or large nicks. Edges should be continuous.

Trials: Articular surfaces should be smooth and free of cracks and deep nicks.

Mating parts: Check to make sure that mating parts fit together without complications. Ensure components are re-assembled securely.

Reamer/drill bits: Inspect "chuck" end for burrs and distortion that might hinder insertion into a drill.

Recommended Sterilization Instructions

Hammering surfaces: Inspect for burrs and large nicks.

Driving devices: Inspect plastic ends for cracks and large nicks.

Metal surfaces: Inspect for corrosion and major deformation.

Powered devices: Verify that power is supplied when the device is turned on and ceases when the device is turned off.

Measuring devices (i.e. gauges, calipers): Check for legible measuring markings.

Maintenance and Care

Prior to sterilization, individual reusable devices that have moving parts (i.e. hinges and sliding parts) may be lubricated using a water soluble lubricant. The use of lubricants is not recommended for implants.

Discard blunt, worn or damaged devices.

Discard devices that do not function properly.

Storage Between Cleaning and Sterilization

Reusable devices that will be stored between cleaning and sterilization should be dry to prevent microbial contamination that could result from wet devices. Containment devices can be stacked for storage. Store in a dry dust-free place.

Recommended Sterilization Instructions

Recommended sterilization methods have been validated to sterility assurance levels (SAL) in compliance with federal and international standards. Other sterilization cycles may also be suitable, but the individuals or hospitals are advised to validate other methods for use with Smith & Nephew orthopaedics products.

Warnings

- It is the responsibility of the user to ensure that the sterilization process is performed using qualified equipment, materials and personnel such that the recommended parameters are achieved.
- Smith & Nephew recommends consulting the Immediate-Use Sterilization statement in AAMI ST 79 for guidance on IUSS.
- Steam is the only method that has been validated for processing by Smith & Nephew. Sterrad or hydrogen peroxide based gas systems have not been validated.
- Package inserts are provided with external fixators to provide directions for resterilization.
- These recommended procedures are intended as a general guide for sterilization of reusable medical devices sold by the Orthopaedics Division of Smith & Nephew. Some devices are labeled with more specific instructions.

- It is the responsibility of the user to validate their sterilization equipment to ensure that the recommended parameters are achieved.

Limitations of Processing

Repeated processing has minimal effect on the devices. Reference the useful life statement from the Recommended Cleaning Instructions on page 6 "Limitations of Reprocessing".

Implants

Implants are sold as single-use devices. The method of initial sterilization for implants that are supplied sterile is noted on the package label. Do not reuse as this may result in product malfunction, failure, or patient injury and may also expose the patient to the risk of infectious diseases.

Containment Devices

Smith & Nephew Orthopaedics containment devices are categorized into families based on the design, density and material of the inner and outer containment devices. Each family has been validated to a 10⁻⁶ Sterility Assurance Level (SAL).

Smith & Nephew designs the proper placement of devices into each containment device layout. The layout has predetermined brackets and artwork descriptions. If devices are added, the user is responsible for validation of the new layout.

When applicable, adhesive labels shipped with the containment device should be applied to the outer containment device.

Large labels typically have a designated outlined box with the instructions "PLACE LABEL HERE." Small labels are placed on the end caps under the labels.

CI and BI placement: Some Smith & Nephew containment device designs have multiple layers of instrument caddies or trays.

Testing has shown that the following CI and BI placement locations should be used for qualification testing by the hospital or healthcare facility:

- **For containment devices that do not contain internal instrument case caddies:** At the center and at each corner of the bottom internal containment device or tray.
- **For containment devices that contain internal instrument case caddies:** At each corner of the caddy on the bottom level of the case.

Recommended Sterilization Instructions *continued*

Containment device weight: Smith & Nephew containment devices are designed to achieve a total case weight (outer case, device, plus inner trays/caddies) of 25 pounds (11.34 kgs) or less. Please be aware that there are some older models that cannot be separated to achieve a weight below 25 pounds (11.34 kgs). These older models have been validated to achieve a 10⁻⁶ SAL.

Only hospital trained personnel should be utilized for inspection and maintenance of containment devices. Modifications to the containment devices should be made only by Smith & Nephew unless the materials and instructions for modifications are supplied by Smith & Nephew (see "Contact Information"). The sterilization of containment devices is validated with the devices placed and positioned in the predetermined placement locations of the containment device. A single absorbent towel (i.e. a huck towel) can be placed under the containment device to aid in drying. Smith & Nephew is the only authorized service/repair company for Smith & Nephew containment devices. Containment devices in need of repair/replacement must be returned to Smith & Nephew. Smith & Nephew does not recommend stacking of containment devices during sterilization.

Preparation for Sterilization

Reusable Medical Devices

It is important that adequate cleaning be carried out prior to sterilization. Reusable medical devices should be placed in suitable packaging for the sterilization process (i.e. central supply wrap [CSR], paper/plastic pouches, rigid containers, etc.) and sterilized prior to each surgical use. Devices are sterilized assembled unless otherwise instructed. If devices were disassembled, ensure components are re-assembled securely with no complications.

Single Use Medical Devices

Trauma plates, nails, screws, pins, and wires are implants and are considered to be single use medical devices. External fixation devices are also considered to be single use medical devices. These devices are sold both nonsterile and sterile and are often removed by the user from their original packaging and placed in a containment device (i.e. instrument case) for processing. These devices should be cleaned prior to sterilization (see the cleaning instructions for Devices **without** challenging design features). These instructions do not apply to processing or reuse of single use devices that have come in contact with blood, tissue or bodily fluids. Prior to sterilization of the device, remove all original packaging and labeling inserts. Place the device in its designated location in the containment device. It is important that adequate cleaning be carried out prior to sterilization.

Sterilization Wrap/Reusable Rigid Containers

Containment devices can be wrapped with an approved CSR wrap or placed in an approved reusable rigid container for sterilization. All sterilization wraps may not be approved for all cycle types. Check with the manufacturer for approvals. Aesculap SterilContainer™ and Case Medical SteriTite® rigid containers with perforated bottoms have been approved for use with Smith & Nephew Orthopaedics instrument sets. These rigid containers are not approved for Immediate Use Steam Sterilization (IUSS) or (Flash Sterilization).

Note to US customers: Sterilizers and wraps used in your sterilization process must be cleared by the FDA.

Water Quality

Use critical water for steam generation.

Recommended Sterilization Instructions *continued*

Recommended Sterilization Parameters

Dynamic Air Removal (Prevacuum) Steam

Exposure temperature:	132° C (270° F)
Exposure time:	4 minutes
or	
Exposure temperature:	135° C (275° F)
Exposure time:	3 minutes
Minimum drying time:	30 minutes

Gravity Displacement Steam

Exposure temperature:	132° C (270° F)
Exposure time:	
Devices not in a containment device	15 minutes
Devices in a containment device	30 minutes*
Minimum vacuum drying:	30 minutes

Immediate Use Steam Sterilization (IUSS)

Exposure temperature:	132° C (270° F)
Exposure time:	4 minutes

*This sterilization cycle is not considered by the United States Food and Drug Administration (US FDA) to be a standard sterilization cycle. Users should only use sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization containers) that have been cleared by the US FDA for the selected sterilization cycle specifications (time and temperature).

For Non-US Customers - UK Steam Cycle Prevacuum cycle

Exposure temperature:	134° C (273° F)
Exposure time:	3 minutes
Vacuum drying:	30 minutes

*World Health Organization (WHO) Steam Cycle***

Exposure temperature:	134° C (273° F)
Exposure time:	18 minutes
Vacuum drying:	30 minutes

**recommended for use in the Russian Federation

Storage After Cleaning and Sterilization

Sterile packaged devices should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, and temperature and humidity extremes. Storage is event related and not time related. Sterile packaged devices can be stored as long as sterile packaging is not breached, or until expiration date. Smith & Nephew does not recommend stacking of wrapped containment devices or rigid containers. Only items sterilized and packaged in materials cleared by the FDA for maintenance of sterility can be stored.

Cleaning and Sterilization Validation Information

Cleaning and Sterilization Validation Information

The detergents, accessories and equipment used to support the validation of cleaning and sterilization were as follows:

Detergents:	Prolystica® 2X Concentrate Enzymatic Presoak and Cleaner, Steris Corporation Klenzyme® Enzymatic Presoak and Cleaner, Steris Corporation
Note: The manufacturer's instructions were followed for use of the detergents.	
Note: Suitability of alternative agents should be checked by the reference to the manufacturer's information and/or physical testing.	
Steam Sterilizer:	AMSCO Model SV-120, AMSCO Lab LV250
Sterilization Accessories:	KC 600 One Step Sterilization Wrap Aesculap SterilContainer™ System Containers JN 440 and JN442 with lid JK489 and Case Medical SteriTite™ SC06FG
Washer-Disinfector:	Steris Reliance 444 Single Chamber Washer/Disinfector Drying Parameters: 116° C minimum for 25 minutes
Ultrasonic Cleaner:	Branson 8510, 40 kHz
Brushes:	Key surgical sterile processing products and O.R. supplier

References

- AAMI TIR 12, "Designing, Testing, and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: a Guide for Device Manufacturers,"
- AAMI TIR 30, "A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices."
- AAMI TIR 34, "Water for reprocessing of medical devices."
- ANSI/AAMI ST79, "Comprehensive guide to steam sterilization and sterility assurance in health care facilities."
- AAMI/ANSI ST81, Sterilization of medical devices- Information to be provided by the manufacturer for the processing of resterilizable medical devices.
- AORN, Perioperative Standards and Recommended Practices.
- Centers for Disease Control, Guideline for Disinfection and Sterilization in Healthcare Facilities, 2008.
- Health Canada June 2011, "Guidance Document on Recognition and Use of Standards under the Medical Devices Regulations."
- ISO 17664, Sterilization of medical devices —Information to be provided by the manufacturer for the processing of resterilizable medical devices.
- German Redbook, Proper Maintenance of Instruments, Instrumenten-Aufbereitung [Instrument Preparation Working Group].
- Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling Guidance for Industry and Food and Drug Administration Staff, March 17, 2015.
- UK Department of Health, Health Technical Memorandum 01-01: Decontamination of reusable medical devices. Parts A and B.
- ANSI/AAMI ST77 "Containment devices for reusable medical device sterilization".
- SO 15883-1, Washers-disinfectors- Part : Requirements for the development, validation and routine control of a sterilization process for medical devices
- ISO 17665-1 Sterilization of healthcare products, most heat- Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices.

Contact Information

Contact Authorized Representative SMITH & NEPHEW, LLC:
Smith & Nephew, LLC, Russia, 105120, Moscow, 2nd
Syromyatnichesky lane., 1, floor 9, prem. 1, room 1 Tel: +7 (495)
755-55-03

Glossary of Symbols

	Medical Device
	Contains hazardous substances
	Sterilized using vaporized hydrogen peroxide
	Single patient – multiple use
 <small>ifu.smith-nephew.com</small>	Electronic instructions of use

Manufacturer	CE Mark	Contact Information
 Smith & Nephew, Inc. 1450 Brooks Road Memphis, TN 38116 USA	 	+1 901 396 2121 Telephone +1 800 821 5700 Information +1 800 238 7538 Orders/Inquiries

European Representative
Smith & Nephew Orthopaedics B.V.
Alemannenstrasse 14
78532 Tuttlingen, Germany
T +1 07462/208-0
F +1 07462/208-135

www.smith-nephew.com

©2015, 2019 Smith & Nephew, Inc.
All rights reserved

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
State of Tennessee)
County of Shelby)

Subscribed (or affirmed) before me this 25th day
of Feb (month), 2021, by Maureen Whitson
(name of signer), whose identity was proven to me by personal
knowledge or satisfactory evidence.

Deborah H. Phillips

(Signature of Notary)

(Seal of Notary)

Total number of pages:

Name: Maureen Whitson MW

Title: Regulatory Affairs Specialist

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