This Apostille is not valid for use anywhere within the United States of America, its territories or possessions.


This Apostille does not certify the content of the document for which it is issued.

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1. Country: United States of America
This public document
2. has been signed by **Karen A. Bushell**
3. acting in the capacity of **Clerk of the Circuit Court for Montgomery County**
4. bears the seal/stamp of the **Circuit Court for Montgomery County**

5. at Annapolis, Maryland
6. the **2nd day of July, 2021**
7. by The Secretary of State of Maryland

9. Seal

John C. Wichman

Secretary of State

State of Maryland, Montgomery County, Sct.

In the Office of the Clerk of the Circuit Court for Montgomery County

I, Karen A. Bushell, Clerk of the Circuit Court for Montgomery County, Maryland, a court of record, hereby certify that **VERA V CHAWLA** was a commissioned/ appointed and qualified Notary Public commencing on the 25th day of August, 2020.

In Testimony Whereof, I have hereunto set my hand and affixed the seal of the Circuit Court for Montgomery County this 29th day of June, 2021.

Karen A. Bushell

Karen A. Bushell

Clerk of the Circuit Court for Montgomery County

2141 P Street NW, Suite 103, Washington, DC 20037

Affidavit

State of Maryland

County of Montgomery

Before me, the undersigned notary public, this day, personally, appeared to me known or proven, who being duly sworn according to law, deposes the following:

I, Yoshiko Yano, hereby certify that the attached document is

- ☒ an original document
☐ a true copy of an original document

to the best of my knowledge and belief.

Signed [Signature]

Subscribed and sworn to before me this

29 day of June, 2024

VERA V. CHAWLA
NOTARY PUBLIC
MONTGOMERY COUNTY
MARYLAND
MY COMMISSION EXPIRES JULY 29, 2024

[Signature]
Notary Public

Info@GlobalRegistrations.US

www.GlobalRegistrations.US

202-909-1660

To whom it may concern,

Date: 6/29/2021

//////////

June 29, 2021

To be submitted to the Competent Health Authorities of the Russian Federation.

1 Bayer Drive
Indianola (PA), United States
15051

www.bayer.com

We, Bayer Medical Care Inc., USA (1 Bayer Drive Indianola, Pennsylvania 15051-0780, USA) hereby state that the attached Instruction for use for MD "XP 65/115VS sterile disposable kit for MEDRAD® MRXperion injection system" is true.

Sincerely yours,

Nishanth Reddy
Sr. Manager, Global Regulatory Strategy Devices - LCE
Bayer Medical Care, Inc.

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Bayer, the Bayer Cross, MEDRAD, MEDRAD MRXperion, MRXperion, MEDRAD FluiDots, and FluiDots may be registered trademarks of Bayer in the US and other countries.

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©2014-2015, 2017, 2019 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.

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A glossary of the symbols used on the MEDRAD® MRXperion sterile disposable products can be found on page 2 of these Instructions for Use.

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ENGLISH

Introduction: Read the information contained in this section. Understanding the information will assist you in operating the device in a safe manner.

Important Safety Notice: This device is intended to be used by individuals with adequate training and experience in diagnostic imaging studies.

Indications for Use: The contents of this package are intended to be used in the delivery of contrast media or saline. They are indicated for single use on one patient only with MEDRAD® MRXperion MR Injectors.

Contraindications: None known.

Restricted Sale: Rx only

Report any serious incident that has occurred in relation to this device to Bayer (radiology.bayer.com/contact) and to your local European competent authority (or, where applicable, to the appropriate regulatory authority of the country in which the incident has occurred).

⚠ Warnings

Re-using fluid containers for more than one procedure may result in biological contamination. Discard contrast and saline containers after filling syringes for a single procedure.

Air embolization can cause death or serious injury to the patient. Do not connect a patient to the injector until all trapped air has been cleared from the syringe and fluid path. Carefully read the instructions for loading and the use of MEDRAD® FluiDots indicators (where applicable) to reduce the chance of air embolism.

Reuse of this product may result in biological contamination, product degradation and/or product performance issues. Properly discard disposable items after single use, or if there is any possibility that contamination may have occurred.

Syringe sterility will be compromised, and patient infection may result, if the plunger is removed from the syringe. Do not remove the plunger to fill the syringe.

Bacterial contamination can occur if syringes are used as storage containers. Use loaded syringes immediately. Do not use syringes to store fluids for later use. Discard unused loaded syringes.

⚠ Warnings

For devices labeled for single use, please note: This product is intended for single use only. Do not resterilize, reprocess or reuse. The disposable devices have been designed and validated for single use only. Re-use of the single use disposable devices pose risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be cleaned or re-sterilized.

Do not use if sterile package is opened or damaged. Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.

Patient or operator injury may result from fluid leaks or syringe or tubing ruptures. Ensure that the fluid path is open and do not use syringes and tubing rated below 350psi (2410 kPa). An occlusion in the fluid path and/or use of syringes or tubing rated below 350psi (2410 kPa) may result in leaks or ruptures.

Patient injury could result if syringe is not properly engaged. Do not load or inject unless the syringe is properly engaged.

Use care in handling and inserting spike into the fluid source. The spike is sharp and may cause personal injury.

Contamination may occur if the end of the spike or luer connector are touched. Do not touch the end of the spike or luer connector.

⚠ Cautions

Component damage or leaks may occur if not installed properly. Ensure all connections are secure; do not over-tighten. This will help minimize leaks, disconnection, and component damage.

Refer to the Injector Operation Manual for further instructions.

Installing New Syringes

NOTE: Use good clinical practice during the loading process to maintain sterility of the syringe kit when opening.

NOTE: Disposable shall be used by use by date as shown on each package.

1. Remove the syringes from packaging.
2. Insert the syringes by pushing them down into the receptacle until they are securely in place (an audible click is heard).
3. Remove the dust covers from the tip of each syringe.
4. Ensure that the plungers are fully advanced.


Figure 1 - 1: Installed Syringes

Filling the Syringes

The syringes on the MRXperion injector can be filled automatically or manually.

NOTE: Use the small syringe (65mL) and small spike to load contrast media into Syringe A. Use the large syringe (115mL) and large spike to load saline into Syringe B.

Filling the Syringes: Automatic

NOTE: During filling, point the injector head up. During an injection, point the injector head down.


NOTE: Auto Fill of Syringe A and Syringe B cannot be started simultaneously.

1. Set and lock protocol with the desired volumes.
2. Install new syringes.
3. Remove the small spike from the packaging. Remove the dust cover from the luer end of the spike. Install the small spike onto the contrast syringe (Syringe A). Do not install with excessive force. Remove the cover from the tip of the spike. Insert the spike into the fluid source.
4. Press the Fill A button twice to begin Auto Fill of Syringe A.

A

5. Ensure the fluid path is free from excess air. Remove the spike.
6. Remove the large spike from the packaging. Remove the dust cover from the luer end of the spike. Install the large spike onto the saline syringe (Syringe B). Do not install with excessive force. Remove the cover from the tip of the spike. Insert the spike into the fluid source.

7. Press the Fill B button twice to begin Auto Fill of Syringe B.

8. Ensure the fluid path is free from excess air. Remove the spike. Discard the fluid source containers and spike.

NOTE: In steps 3 and 6, if using the Female/Female Adapter (FFA) to load fluid from a pre-filled syringe, attach the FFA onto the tip of the syringe instead of a spike. (The FFA may be purchased separately. The catalog number for this item is FFA 50.)

9. Connect the disposable tubing following the instructions in the **Connector Tube Installation** section.

Filling the Syringes: Manual

NOTE: During loading, point the injector head up. During an injection, point the injector head down.

1. Install new syringes.
2. Remove the small spike from the packaging. Remove the dust cover from the luer end of the spike. Install the small spike onto the contrast syringe (Syringe A). Do not install with excessive force. Remove the cover from the tip of the spike. Insert the spike into the fluid source.
3. Fill Syringe A by pressing the Enable Piston Control button to enable the piston controls and use the reverse piston control to fill the syringe with the desired amount of fluid. Alternatively, use the manual knob on Side A of the injector.
4. Ensure the fluid path is free from excess air. Remove the spike.
5. Remove the large spike from the packaging. Remove the dust cover from the luer end of the spike. Install the large spike onto the saline syringe (Syringe B). Do not install with excessive force. Remove the cover from the tip of the spike. Insert the spike into the fluid source.
6. Fill Syringe B by pressing the Enable Piston Control button to enable the piston controls and use the reverse piston control to fill the syringe with the


desired amount of fluid. Alternatively, use the manual knob on Side B of the injector.

7. Ensure the fluid path is free from excess air. Remove the spike. Discard the fluid source containers and spike.

NOTE: In steps 2 and 5, if using the Female/Female Adapter (FFA) to load fluid from a pre-filled syringe, attach the FFA onto the tip of the syringe instead of a spike. (The FFA may be purchased separately. The catalog number for this item is FFA 50.)

8. Connect the disposable tubing following the instructions in the **Connector Tube Installation** section.

To help avoid an air embolism, syringes from Bayer are equipped with FluiDots indicators. FluiDots indicators should be observed as part of an arming procedure. When the FluiDots are viewed through an empty syringe, the dots appear as small narrow ellipses. When viewed through a fluid-filled syringe, the dots become larger (almost round).


Operator vigilance and care, coupled with a set procedure, is essential to minimizing the possibility of air embolism. Ensure that one operator is designated the responsibility of filling the syringe(s). Do not change operators during the procedure. If an operator change must occur, ensure that the new operator verifies that the fluid path is purged of air.

Connector Tube Installation

After loading and priming the syringes, install the connector tube.

1. Remove the connector tube from the package. Remove the dust covers on the luer fittings.
2. Ensure that all air is purged from the syringes.
3. Securely attach the connector tube to the syringes as shown. Do not install with excessive force.
4. Ensure that the connector luer fitting is secured to the tips of the syringes and verify that the tubing is not kinked or obstructed.
5. Prime the tubing by pressing the Prime button.
6. Ensure that all air is purged.

7. Rotate the injector head downward.

8. Examine the fluid path for air. Confirm that the fluid path is clear of excess air by pressing the Check for Air Confirmation Button. Connect to the patient.

Removing the Syringes

1. Disconnect the disposable tubing set from the vascular entry device. The disposable tubing set does not need to be disconnected from the syringes.
2. Rotate the syringes and gently pull the syringes out of the injector head. Discard the syringes and disposable tubing set.

NOTE: Once the syringes are removed from the injector head, the piston will automatically retract.

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Certified

5. at Annapolis, Maryland
6. the **2nd day of July, 2021**
7. by The Secretary of State of Maryland
8. **No. 560409**
9. Seal

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I, Karen A. Bushell, Clerk of the Circuit Court for Montgomery County, Maryland, a court of record, hereby certify that VERA V CHAWLA was a commissioned/ appointed and qualified Notary Public commencing on the 25th day of August, 2020.

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

1 Bayer Drive
Indianola (PA), United States
15051

www.bayer.com

We, Bayer Medical Care Inc., USA (1 Bayer Drive Indianola, Pennsylvania 15051-0780, USA) hereby state that the attached addendum to the Instruction for use for MD "XP 65/115VS sterile disposable kit for MEDRAD® MRXperion injection system" is true.

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Addendum to the instruction for use

1. Operation conditions

	Conditions
Temperature	10°C to 35°C
Humidity	20% to 90%, non-condensing
Ambient pressure	from 70 kPa to 103 kPa

2. Transportation and storage conditions

Product	Conditions
Temperature	-20°C to 60°C
Humidity	5% to 100%
Ambient pressure	from 57 kPa to 106 kPa

3. Technical characteristics

Parameter	Value
Syringe for contrast media, 65 ml	
Volume, ml	65
Length, mm	196.5
Syringe for saline solution, 115 ml	
Volume, ml	115
Length, mm	200.6
Small needle for contrast media	
Needle length without cap, mm	35,81
Length of the protective cap, mm	34,8
The large needle for saline solution	
Needle length, mm	68,07
Low pressure connecting tube with T-shaped connector	
Straight tube length, mm	- 190,5
Twisted tube length, mm	2489,2
Straight tube outer diameter, mm	2,54
Twisted tube outer diameter, mm	3.18

Page 3 of 3

4. Shelf life

Sterilization of set is performed by gamma radiation.

Shelf life (sterilization validity): 3 years from the date of manufacture/sterilization.

5. Disposal Procedure

Used kit components must be disposed of in a biohazard waste container in accordance with common medical practice and local regulations.

6. Warranty liabilities

The warranty lifetime is not provided for due to the specifics of the medical device.

Manufacturer warrants that the products are free from defects for the entire device shelf-life.

Manufacturer bear no responsibility nor liability in case of device misuse or non-adherence to the Instructions for use.

7. Reclamation

Issues related to the medical device quality must be referred to the manufacturer's representative in the Russian Federation at the address below:

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