
Legal Manufacturer:

Siemens Healthcare GmbH
Henkestr. 127, 91052 Erlangen
Germany

User manual Supplement for Russia
from April 2022

***The magnetic-resonance system MAGNETOM Aera
with accessories***

This document has been printed by the manufacturer in Erlangen, Germany.

Friederike Hertle
Head of Regulatory Affairs MR
Siemens Healthcare GmbH, Magnetic Resonance (MR)

1. The name of the Medical Device: the accurate of the MD.

The magnetic-resonance system MAGNETOM Aera with accessories.

2. Information about Manufacturer of the Medical Device:

2.1 Legal Manufacturer

Siemens Healthcare GmbH
Address: Henkestr. 127, 91052 Erlangen, Germany
Telephone: +49 9131 84-0
Web: www.siemens-healthineers.com

2.2 Manufacturing site

Siemens Healthcare GmbH, Magnetic Resonance (MR)
Address: Henkestr. 127, 91052 Erlangen, Germany
Telephone: +49 9131 84-0
Web: www.siemens-healthineers.com

2.3 The name and address of the local representative in Russia

Siemens Healthcare LLC
Address: Dubininskaya Str. 96, 115093 Moscow, Russia
Tel.: +7 495 737 12 52
Web: www.siemens-healthineers.com/ru

3. Designation of the medical device: Purposes for use (Medical Device operation data):

3.1 Intended Use

The MAGNETOM Aera system is indicated for use as a magnetic resonance diagnostic device (MRDD) that produces transverse, sagittal, coronal and oblique cross sectional images, spectroscopic images and/or spectra, and that displays the internal structure and/or function of the head, body, or extremities.

Other physical parameters derived from the images and/or spectra may also be produced. Depending on the region of interest, contrast agents may be used. These images and/or spectra and the physical parameters derived from the images and/or spectra when interpreted by a trained physician, yield information that may assist in diagnosis.

The MAGNETOM Aera system may also be used for imaging during interventional procedures when performed with MR compatible devices such as in-room display and MR Safe biopsy needles.

3.2 Indications

Magnetic Resonance Imaging (MRI) is a noninvasive technique used for diagnostic imaging. MRI with its soft tissue contrast capability enables the trained user to differentiate between various soft tissues, for example, fat, water, and muscle, but can also discriminate between the different bone structures. This allows for visualization of anatomic details and can offer more sensitivity and specificity in the detection of certain abnormalities than other imaging modalities, when read and interpreted by trained personnel. In general, it is the health care provider who decides on the appropriate diagnostic procedure based on, for example, the patient's physical examination and medical history.

3.3 Contraindications

In general, an MR examination is contraindicated for patients with electronic or electronically conductive implants or metals, especially those containing ferromagnetic foreign matter.

Typical contraindications for MR examinations are:

- Electronic implants: e.g. pacemakers, stimulators, insulin pumps
- Artificial heart valves, aneurysm clips
- Metal splinters in the eye (danger of retinal detachment)
- Artificial anus (anus praeter) with magnetic closure
- Transdermal drug patches with metallic backings
- Electrically conductive implants and prostheses
- Metallic spirals for contraception (IUDs = Intra Uterin Devices)
- Transdermal or other similar implants for body piercings as well as magnetic piercings
-

3.4 Possible side effects

Possible undesirable side effects for MR are dizziness, heating, claustrophobia and nerve stimulation.

3.5 Mode of application (Medical Device application data)

1. Preparing an examination
 - Patient registration
 - Positioning the patient on the patient table
 - Positioning the local coils and cushions
 - Table positioning
2. Performing an examination
 - Measuring of the Localizer and parallel loading of previous reference images
 - Planning the examination volume
 - Scanning
 - Viewing and editing images during examinations
3. Basic syngo MR viewing, post-processing, and filming functionalities
 - Viewing
 - Mean Curve
 - Composing/Combine
 - 3D
 - Filming

3.6 Conditions for use

For ambient conditions in hospitals or MRI medical practices please refer as following:

Control room

Temperature	15°C–30°C (59°F–86°F)
Relative air humidity	40%–60%
Absolute air humidity	< 11.0 g/kg

Control room

Air pressure	1060 hPa–700 hPa ²⁴⁾
--------------	---------------------------------

Equipment room

Temperature	15°C–30°C (59°F–86°F)
Relative air humidity	40%–80%
Absolute air humidity	< 11.0 g/kg
Air pressure	1060 hPa–700 hPa ²⁵⁾

Examination room

Temperature	18°C–22°C (64°F–72°F)
Relative air humidity	40%–60%
Absolute air humidity	< 11.0 g/kg
Air pressure	1060 hPa–700 hPa ²⁶⁾

4. Classification of the Medical Device:

4.1 Type of the MD in accordance with the nomenclature for classification

• **Europe (CE)**

Kind of product:	Active medical device for diagnosis
Classification rules according to MDD 93/42/EEC, Annex IX:	Rule 10
Classification according to Annex IX:	Ila
Classification Criteria:	Active devices intended for diagnosis are in Class Ila if they are intended to supply energy which will be absorbed by the human body, except for devices used to illuminate the patient's body, in the visible spectrum

• **Russia**

Classification in accordance with Russia Nomenclature Classification	135160
Class A - Waste generated during / after the use of medical devices must be classified and disposed of (disposed of) as epidemiologically safe waste, close to MSW in composition, in accordance with current legislation.	

4.2 Software classification according to IEC 62304:2006.

Based on the current product risk management results, the SW safety class is "B". This classification is subject to review whenever changes in this risk analysis document apply.

4.3 Class of the potential risk of the medical device application

Potential risk in accordance to SW safety class "B": non-serious injury possible.

Brand name and version of SW:

syngo MR E11

syngo MR XA30A

4.4 Type of contact with the human body

MR system device is not intended to be sterile.

4.5 Laser classification

The laser of the laser light localizer is classified as Class 1M according to IEC 60825-1 (Class II according to US CDRH).

4.6 IP Protection class

IP protection class according to IEC 60529 is IP X0.

4.7 Applied parts classification

Patient table	Type B
Local coils	Type B or BF
Body coil	Type B
ECG/Pulse module	Type BF

4.8 Operating mode

Continuous operation with short-term load.


5. Technical Specifications of the Medical Device (the main technical parameters of the Medical Device; the MD's scheme (design), characteristics, properties, features):

5.1 The photo of the device

Below you will find the front view and side view picture of the MAGNETOM Aera.


Front view with fixed patient table - photo

Front view with fixed patient table –scheme in mm


Side view with fixed patient table - photo


Side view with fixed patient table –scheme in mm

5.2 The list of components of the medical device


No.	Name of component	Description	Photograph
I. Basic block:			
1.	MAGNETOM Aera – System	Whole-body superconductive Zero Helium Boil-Off 1.5T magnet.	Standard variant: Pure White Design with brilliant front ring 
2.	Gradient system Tim [204x24] XJ Gradients Tim [204x48] XJ Gradients Tim [204x48] XQ Gradients Tim [204x64] XQ Gradients	The whole-body gradient system engine.	
			 GPA 62/20/F (XJ gradients) GPA 90/22/F(XQ gradients)
			
			

No.	Name of component	Description	Photograph
I.	Basic block:		
3.	Transmit/receive Body Coil	Integrated no tune transmit/receive Body Coil	
4.	Magnet cover	Magnet cover	Standard variant: Pure White Design with brilliant front ring  Illumination MoodLight™ variant with brilliant front ring  Water Blue Design with brilliant front ring


No.	Name of component	Description	Photograph
I.	Basic block:		
			 Light Green Design with brilliant front ring with customized Logo 
5.	Integrated magnet cooling system	Integrated magnet cooling system	Chiller 45kW 098  Chiller 60kW 098  Interface Panel IFP 098



No.	Name of component	Description	Photograph
I.	Basic block:		
6.	Patient table	The patient table is mounted directly to the magnet assembly. The table comes with a comprehensive set of patient positioning cushions. The Combi Dockable Table is specially designed for flexible patient positioning. It only requires repositioning once on a transfer plate, transporting the patient safely to the MR system.	 Fix table 1.5T  Mobile table 098 1.5T  Combi Dockable Table
7.	Electronics cabinets	The electronics cabinet contains the electronics for Measurement sequence control generating the magnetic field gradients radio frequency generation Voltage supply for network Integrated cooling system The system shall support two cooling scenarios: passive and active cooling.	 (1) Gradient cabinet (2) Control cabinet (3) System separator


No.	Name of component	Description	Photograph
I.	Basic block:		
8.	Alarm box	The alarm box has the following functions: - Displays alarm signals - Switches the MR system on and off - Magnet Stop / Magnet Quench The alarm box is installed near the syngo Acquisition Workplace.	
9.	Emergency switches	The MR system has different types of emergency switches. - Magnet Stop Shutting down the static magnetic field (quenching) - in case of accidents with attracted metal parts or in case of fire - Emergency Shutdown Electric power of the entire MR system is switched off, but magnet remains at field, e.g. in case of fire - Table Stop Motorized table movement is stopped. In case of accidents or injury due to table movement	Magnet Stop  Emergency Shut-down switch  Table Stop Button On the intercom  On the patient table 


No.	Name of component	Description	Photograph
I. Basic block:			
10.	Intercom	The intercom allows personnel and patients to communicate during the examination. In addition, some important operations like stopping the patient table can be managed from the intercom. Optionally, music or automatic voice outputs can be played in the examination room via the loudspeaker or the headphones.	 (1) Stop button (2) Talk button (3) Reset Patient Alert button (4) Listen button (5) Physio Signal button (6) Audio connection (7) Reset Stop button
11.	Acquisition Workplace	The syngo Acquisition Workplace (syngo Acq WP). It includes the host processor with the operating elements monitor, keyboard, and mouse.	
12.	Keyboard	The syngo Acquisition Workplace is equipped with an original Siemens keyboard. This keyboard is a modified Windows keyboard where the numeric keys have been replaced with symbol keys.	
13.	Head/Neck Coils Head/Neck 16 MR Coil	The Head/Neck coil supports high quality imaging. The unique DirectConnect	Head/Neck 16 MR Coil


No.	Name of component	Description	Photograph
I.	Basic block:		
	Head/Neck 20 MR Coil	technology allows users connecting elements of the Head/Neck coil without cables.	 Head/Neck 20 MR Coil 
14.	Spine Coils Spine 32 MR Coil Spine 24 MR Coil	The Spine 32 is used for examining the spine and for combined measurements with other coils. The 24-channel coil with its 24 integrated pre-amplifiers ensures excellent signal-to-noise ratio. The unique DirectConnect technology allows connecting the 24 coil elements of the Spine 24 without the need to plug in any cable.	Spine 32 MR Coil  Spine 24 MR Coil 
15.	Body Coils Body 18 MR Coil Body 18 long MR Coil Body 6 MR Coil Body 6 long MR Coil	The body coil with its integrated pre-amplifiers ensures excellent signal-to-noise ratio. The SlideConnect plug allows for fast and easy patient preparation resulting in less table time. Fast acquisition times enabled by iPAT in all directions.	Body 18 MR Coil / Body 18 long MR Coil 

No.	Name of component	Description	Photograph
I.	Basic block:		
		Both coils are optionally available with long coil cable.	Body 6 MR Coil / Body 6 long MR Coil 
16.	Flex Coil Interface	The Flex Coil Interface is required for connecting the Flex coils, the Loop coils, the Special Purpose coil and the Endorectal coil.	
17.	Flex Large 4 Coil	Flex Large 4 Light-weighted, very flexible, iPAT compatible, 4-element no-tune receiver coils which are made of soft and smooth material.	Flex Large 4 MR coil 
	UltraFlex Large 18 Coil	The coil flexible large (UltraFlex Large 18) is used for inspection of large extremities (for example, from an average to a big shoulder joint, a basin, a knee joint, an ankle and hands), and also for inspection of an abdominal cavity.	UltraFlex Large 18 coil 
18.	Flex Small 4 Coil	Flex Small 4 Light-weighted, very flexible, iPAT compatible, 4-element no-tune receiver coils which are made of soft and smooth material.	Flex Small 4 MR coil 
	UltraFlex Small 18 Coil	The coil flexible small (UltraFlex Small 18) is	UltraFlex Small 18 coil



No.	Name of component	Description	Photograph
I.	Basic block:		
		used for inspection of small extremities (for example, from small to an average shoulder joint, a small ankle, an elbow and hands), and also for inspection of an abdominal cavity.	
19.	Tx/Rx CP Head Coil for vascular examinations	Tx/Rx CP Head Coil is suitable for the following examinations: - Vascular examinations in the head - e. g., vascular examinations of the base of the skull and the cerebrum	Tx/Rx CP Head Coil 
20.	Body Coils for expert examinations Body 30/60 MR Coil Contour 24/48 Coil	Flexible, universal receive coil with integrated preamplifiers	Body 30 MR Coil  Body 60 MR Coil (The Body 60 coil represents a combination of two Body 30 coils)  Contour 24 coil 

No.	Name of component	Description	Photograph
I.	Basic block:		
			 Contour 48 coil (The Contour 48 coil represents a combination of two Contour 24 coils) 
21.	Peripheral Angio 36 Coil	The coil for a peripheral angiography (Peripheral Angio 36) is suitable for inspections of vessels of the lower extremities. Depending on requirements the coil can be used in provision by the head forward or legs forward. The Peripheral Angio 36 coil can be used together with the coil for a backbone and one or more coils for a body. Thanks to a combination of these coils the field of inspection extends.	Peripheral Angio 36 Coil 
22.	Tx/Rx Knee 15 Coils for joint examinations	Transmitter/receiver coil for joint examinations in the area of the lower extremities.	Tx/Rx Knee15 Flare MR Coil 

No.	Name of component	Description	Photograph
I.	Basic block:		
			Tx/Rx Knee 15 MR Coil 
23.	Tx/Rx Knee 18 Coil	The coil for a knee joint (TxRx Knee 18) is used for examinations of the left or right knee. The coil can be moved on a basis plate on a longitudinal axis to bring her in the center for receiving a picture of optimum quality.	Tx/Rx Knee 18 Coil 
24.	Tx/Rx CP Extremity Coil	The TxRx CP Extremity is used for examining knee joints and foot ankles. This coil is available for 1.5T MR systems only. The coil can be shifted laterally on the base plate to center it for optimal image quality. There are two predefined positions on the patient table to position the base plate (knee position and foot ankle position).	Tx/Rx CP Extremity Coil 
25.	Foot/Ankle 16 Coil	The Foot/Ankle 16 is used to examine feet and ankles	Foot/Ankle 16 Coil 

No.	Name of component	Description	Photograph
I. Basic block:			
26.	Hand/Wrist 16 Coil	The Hand/Wrist 16 is used to examine wrists, hands, and fingers. Depending on the requirements, the coil may be used with the patient in the head first (prone position) or feet first (supine position) orientation. For the feet first orientation, the coil is positioned on the base plate.	Hand/Wrist 16 Coil 
27.	Shoulder Large 16 Coil Shoulder Small 16 Coil	The Shoulder Small/Large is used for examining the shoulder. The coil is used in the head first orientation.	Shoulder Large/Small 16 Coil 
28.	Shoulder Shape 16 Coil	The coil for expert examinations of a shoulder joint (Shoulder Shape 16) is used for inspection of a shoulder and is supplied with a flexible shoulder cup.	Shoulder Shape 16 Coil 
29.	Breast 18 coil for expert examinations	The Breast coil is used for examining female patients. Depending on the requirements of the examination, the coil may be used in the head first or feet first orientation. Both sides of the breast are displayed simultaneously to allow for comparative evaluation.	Breast 18 Coil 



No.	Name of component	Description	Photograph
I.	Basic block:		
30.	Breast Biopsy Coil 4-channel BI Breast Coil	The coil can be used for biopsies or high-resolution breast imaging. It consists of two parts that are positioned in the frame.	4-ch BI Breast Coil 
	Breast Coil for visualization and biopsy Breast BI 7 Coil	The Breast BI 7 coil can be used for visualization and a biopsy of a mammary gland. Depending on the requirements of the examination, the coil can be used in head first or feet first orientation.	Breast BI 7 Coil 
31.	Modular Breast Coils for visualization and biopsy 2-/4-/8-ch Sentinelle Breast Coil 2-/10-/16-ch Sentinelle Breast Coil	The Sentinelle Breast Coil supports the Grid biopsy method. The Sentinelle Breast Coil delivers brilliant image quality for high-resolution 2D and 3D MR breast imaging. Techniques for reducing scan times, such as parallel imaging, can be used very well.	2-/4-/8-Channel Sentinelle Breast Coil  2-/10-/16-Channel Sentinelle Breast Coil 




No.	Name of component	Description	Photograph
I. Basic block:			
32.	Prostate 2 Coil	The Prostate 2 is used for examinations of the prostate, colon and the surrounding tissue in the pelvis. Depending on the requirements of the examination, the coil may be used in the Head First or Feet First orientation.	Prostate 2 Coil  (1) Prostate 2 coil (2) Mounting Clip (3) Stabilization Arm (4) Patient support
33.	Special Purpose 4 Coil	The special purpose coil is a no-tune receiver coil designed for examination with small Field-of-Views.	Special Purpose 4 Coil 
34.	Loop Coil (small)	No-tune loop receive coil for examinations of small structures near the surface (e.g. hands, fingers, skin).	Loop 4 Coil (small) 
35.	Loop Coil (medium)	No-tune loop receive coil for examinations of inner ear, structure of wrist and fingers, pediatric examinations.	Loop 7 Coil (medium) 
36.	Loop Coil (large)	No-tune loop receive coil for examinations of the upper or lower extremities (e.g. shoulder).	Loop 11 Coil (large) 
37.	Prostate Endo Coil for 1.5T MR	The Medrad disposable Prostate coil, through its high-resolution imaging and small field of view, provides improved ability to visualize internal architecture of the prostate and periprostatic structures, including prostate capsule and vascular bundles. The radiological images may be used in	Prostate Endo Coil 


No.	Name of component	Description	Photograph
I.	Basic block:		
		diagnosing diseases of the prostate. The inflated balloon makes a positioning of the Prostate coil in close proximity to the wall adjacent to the prostate gland possible.	
3.	Colorectal Endo Coil for 1.5T MR	The Medrad disposable Colorectal coil, through its high-resolution imaging and small field of view, provides images of the colorectal anatomy, which may be used to aid in diagnosing colorectal diseases.	Colorectal Endo Coil 
39.	Cervix Endo Coil for 1.5T MR	The Medrad disposable Cervix coil, through its high-resolution imaging and small field of view, provides images of cervix, uterine, and endometrial anatomies, which may be used to aid in diagnosing diseases in these anatomical regions.	Cervix Endo Coil 
40.	Pediatric 16 Coil	16-channel receive coil for head and neck imaging of newborns and children up to 18 months of age.	Pediatric 16 Coil 




No.	Name of component	Description	Photograph
I. Basic block:			
41.	Workstation-based Solution	syngo.via Workstation is a syngo.via-based entry level sales package for advanced visualization customers. It provides the same software functionality, as the conventional syngo.via product.	
42.	Dot Control Units	The ergonomically designed Dot Control Units are integrated left and right into the front covers for controlling table movement and interaction with the Dot Display.	
43.	Connection cables	Connection cables	
1.	Measurement phantoms	Measurement phantoms are used for quality measurements.	Plastic bottle, 15ml  Plastic bottle, 1900ml  Plastic bottle, 5300ml 


No.	Name of component	Description	Photograph
I.	Basic block:		
			Spherical phantom D165  Spherical phantom D240  Hand/Wrist 16 phantom  Foot/Ankle 16 phantom  Special phantom for Endorectal 
45.	Phantom holders	Phantom holders	Phantom holder 17 


No.	Name of component	Description	Photograph
I.	Basic block:		
			Phantom holder 18  Phantom holder 23  Phantom holder 24  Phantom holder 25  Phantom holder 26  Phantom holder 27  Phantom holder 28 

No.	Name of component	Description	Photograph
I.	Basic block:		
			Phantom holder 29  Phantom holder 31  Phantom holder 32  Phantom holder 40  Phantom holder 45  Phantom holder 46  Phantom holder 47 



No.	Name of component	Description	Photograph
I.	Basic block:		
			Endorectal coil simulator  Prostate 2 phantom holder 
46.	The application package "Tim Application Suite": - Basic software license: basic license key for running the syngo software; - Basic MR applications license; - 3D functionality; - Argus (cardiac) basic viewing; - DICOM worklist support; - PhoenixZIP; - Different standard measurement/examination and post-processing functionality.	The Tim Application Suite offers a complete range of clinically optimized sequences, protocols and workflow functionalities for all body regions.	Neuro Suite  Angio Suite  Cardiac Suite  Body Suite 

No.	Name of component	Description	Photograph
I.	Basic block:		
			 Onco Suite  Ortho Suite  Breast Suite
47.	MAGNETOM Aera – System Russian User Manual	Instruction for use.	 
48.	Physiological Measurement Unit	Physiological Measurement Unit (PMU) - Wireless Physio Control for wireless triggering, synchronizes the measurement with the physiological cycles of cardiac and/or respiratory motion. Additional PERU, PPU and charging stand.	

No.	Name of component	Description	Photograph
I.	Basic block:		
49.	MAGNETOM Family – OM MR system	Instruction for use.	
50.	SOM – MAGNETOM Aera	Instruction for use.	
51.	MAGNETOM Family – OM - Coils	Instruction for use.	
52.	OM - scanning and post-processing	Instruction for use.	

No.	Name of component	Description	Photograph
I.	Basic block:		
53.	MAGNETOM Aera – OM Implant Suite	Instruction for use.	
54.	MAGNETOM Family – OM - 2-/4-/8-Channel Sentinelle Breast Coil, 2-/10-/16-Channel Sentinelle Breast Coil	Instruction for use.	
55.	MAGNETOM Family – OM - 4 Ch BI Breast / 16 Ch AI Breast	Instruction for use.	
56.	MAGNETOM Family – OM Breast BI 7	Instruction for use.	

No.	Name of component	Description	Photograph
I.	Basic block:		
57.	MAGNETOM Family – OM Contour 24 Coil	Instruction for use.	
58.	MAGNETOM Artis Combi Suite – Combi Dockable Table	Instruction for use.	
59.	MAGNETOM Combi Suite Neurosurgery – Combi Dockable Table (Maquet Hybrid)	Instruction for use.	
60.	MAGNETOM Combi Suite Neurosurgery – Combi Dockable Table (Trumpf Medical)	Instruction for use.	

No.	Name of component	Description	Photograph
I.	Basic block:		
61.	SOM - System Maintenance Schedule MAGNETOM Aera/Aera XJ/Skyra/Avanto fit/Skyra fit	Instruction for use.	
62.	MAGNETOM Family Addendum Regulatory Information	Instruction for use.	
63.	Transfer board	The transfer board for neurosurgical applications with head holders for intraoperative MRI	Transfer board Maquet  Transfer board Trumpf  Transfer board Artis 



No.	Name of component	Description	Photograph
I. Basic block:			
64.	RT Positioning Package CIVCO 3002	Package includes: 50cm wide Overlay, Lok-Bar, 2 each Body Bridges and 1 each Head Bridge.	
5.	NORAS Headholder CombiSuite	Starter package for interventional neurosurgical MRI, consisting of a NORAS head holder in combination with a NORAS OR head coil system, consisting of a 4ch top array and a 4ch bottom array.	
66.	Patient Transportation devices	MRI-Stretcher, non magnetic, chrome-nickel compounded for patient transport inside the clinic building. Continuously adjustable head piece with spring brake, synthetic leather cushion, movable barrier on either side, stretcher movable on four fixable wheels diameter 100 mm (for max. 3 mm floor unevenness). Dimensions: 190 x 65 x 80 cm (l x w x standard height) Maximum laminary load: 200 kgWeight (empty): 41 kg Standard upholstery colour: blue Additionally available upholstery colours:	



No.	Name of component	Description	Photograph
I.	Basic block:		
		Red, green, birch, grey, yellow, white, violet, ice, light brown, light blue, coffee brown, turquoise blue, black, light grey, magenta, dark-pink, tangerine, atoll blue, rosewood and charcoal-grey.	
7.	IV Stand, non-magnetic	IV-Rod, fix mounted, for AGA non magnetic stretcher. - Diameter: 180mm - Length: 1720mm - With four hooks	
68.	ECG-electrods	Disposable electrodes for long-term monitoring - For children and adults - Neoderm support - Gelled with solidgel - Carbon coated - X-ray translucent - Square electrodes with rounded edges - Diameter: 51mm	

No.	Name of component	Description	Photograph
II.	Accessories:		
1.	syngo MR Workplace	The syngo MR Workplace allows for evaluation, documentation and post-processing of previously measured images while acquiring images at the syngo Acquisition Workplace.	
	Keyboards	An original Siemens keyboard. This keyboard is a modified Windows keyboard where the numeric keys have been replaced with symbol keys.	
3.	Foot switch In-Room	MR Foot Switch In-Room with two pneumatic buttons for start and stop of a preset MR sequence. The foot switch is MR compatible and is positioned near the patient table on the examination room floor. The MR foot switch is especially suited for sequence control when MR examinations are carried out during operations.	
4.	Dual Monitor Package	The Dual Monitor Option will improve the ease of use of MR examination, with a more smooth and comfortable work space that avoids interruptions between planning, scanning and viewing. It allows to keep running patient examinations always	

No.	Name of component	Description	Photograph
II.	Accessories:		
		in sight to allow for fast interactions. Additional to the host monitor a second monitor of the identical type is included in this option, being connected to the second port of the host graphics card. For all tasks on both monitors still one mouse and keyboard be used. The left monitor will be used for the examination exclusively, whereas all postprocessing tasks are performed on the right monitor.	
5.	Gradient supervision	The supervision monitors that cables, connections, or other components of the gradient system do not show excessive heating.	
6.	Application cushion	Application cushion	
7.	Respiratory cushion	Respiratory cushion with pressure hose	

No.	Name of component	Description	Photograph
II.	Accessories:		
8.	Respiratory belt	Respiratory belt	
9.	Standard Headphone	The standard Siemens headphones are intended for communication with the patient and can be used in combination with ear plugs.	
10.	The MagnaCoil Headset (In-Ear Headphone)	The MagnaCoil Headset can be used as an alternative headphone (for example, in head coils) to the standard Siemens headphones which are part of the products since years. Sometimes the standard headphones are not very practicable and therefore Siemens wants to offer an alternative: the MagnaCoil Headset. This In-Ear Headphone works for easy communication with patient while using the head coil. The user can communicate better with the patient during the exam. The patient can listen to music	
11.	In-Room syngo Acquisition Workplace - an additional	The In-Room syngo Acquisition Workplace is an additional operating console in	





No.	Name of component	Description	Photograph
II.	Accessories:		
	operating console in the examination room	the examination room. It is used for image viewing and MR system operation and comprises the following components: - LCD monitor - Tray with trackball and keys - Foot switch (optional)	
12.	Positioning frame for Breast Biopsy Coil 4-ch BI Breast Coil	Open design positioning frame for inserting the 4-channel BI Breast Coil. The positioning frame comes supplied with: - a head cushion - fixation rail to proper position and fix the coil positioning frame on the patient table - abdominal support cushions - Breast stabilizing cushions - positioning guidance	
13.	Breast Biopsy Starter Kit for Breast Biopsy Coil 4-ch BI Breast Coil	The breast biopsy kit contains a fixation unit plus accessories for both the Post-Pillar and the Grid biopsy method. - Insertion plate; - Horizontal and vertical slat plates; - Needle guide post/pillar method - Height adjustable grids reusable and/or disposable - Frame fixture This is a training kit, not used with patients.	Insertion plate  Horizontal and vertical slat plates  Needle guide (post/pillar method) 









No.	Name of component	Description	Photograph
II.	Accessories:		
			 Height adjustable grids reusable and/or disposable  Frame fixture 
14.	Coil adapter	Coil adapter plug for up to 8 receive and 1 transmit channels, in order to connect existing dedicated knee and breast coils (Tx/Rx 15-channel Knee Coil, CP Extremity Coil, 4-channel BI Breast Coil, 16-channel AI Breast Coil, (2/4)/8-channel Sentinelle BreastCoil and (2/10)/16-channel Sentinelle BreastCoil) with all MAGNETOM 1.5T Systems using Tim 4G-technology	 (1) Socket cover for coil socket (2) Release for coils (3) Interface coil plug (4) Release for Tim Coil Interface
15.	Multi-channel Interface	The Multi-Channel Interface transforms the Head Coil Direct Connect sockets into two Tim 4G Slide Connect sockets. In	


No.	Name of component	Description	Photograph
II.	Accessories:		
		combination with Tim Coil Interfaces (not included in this option), this allows to connect specific Legacy or 3rd party coils.	
16.	fMRI Trigger Converter	An optical trigger signal is available to trigger external stimulation devices in fMRI experiments.	
17.	Patient Supervision TV #T+D Additional Camera #T+D Supervision Camera & Controller #T+D	This package contains a special video camera for monitoring the patient during an MR examination. The supervision solution is customizable and designed to address different site specific requirements. Up to 4 cameras can be optionally connected for patient supervision in the examination or waiting room. Additional video camera to allow for comprehensive monitoring of the patient's face independent of patient position during an MR examination. Wall mount for the patient monitor.	






No.	Name of component	Description	Photograph
II.	Accessories:		
18.	Patient TV wall support Patient TV wall support #T+D	Wall mount for the patient monitor	
19.	LCD Color Remote Monitor #T+D	High-resolution LCD color flat screen monitor, with integrated gamma correction for optimum display of MR grayscale images and automatic backlight control for long-term brightness stability.	
20.	INSIGHTEC HIFU Support	The package supports customers to perform MR-guided HIFU procedures with the INSIGHTEC Exablate system. It includes the following hardware components: - hardware router - trigger converter as well as 3 software licenses: - Access-i - Phase Imaging - Projection	Router  Trigger Converter 






No.	Name of component	Description	Photograph
II.	Accessories:		
21.	Dot display	The Dot (Day Optimized Throughput) display provides status information as well as several tools. The display is located above the magnet bore opening at the front of the magnet cover. Optionally, at the back of the magnet, a second display can be installed for the additional control unit.	
22.	Charging station	The charging station is installed separately near the syngo Acquisition Workplace and is used for storing both units.	
	Additional PMU Sensor Kit	Additional PERU, PPU and charging stand.	Sensor Charging Station 098  PERU_098 Phys. ECG/Respiratory Unit PPU_098 Peripheral Pulse Unit





No.	Name of component	Description	Photograph
II.	Accessories:		
			
24.	Cushions and positioning aids	Cushions and positioning aids - Table pad 1 - Table pad 2 - Table pad 3 - Table pad 4 - 4 belt straps - 2 different-sized sandbags, 2 each - Head-heel insert - Lookout mirror holder	 Table pad 1  Table pad 2  Table pad 3  Table pad 4  Head-heel insert  4 belt straps  2 different-sized sandbags, 2 each  Lookout mirror holder
25.	Set of positioning aids	Set of positioning aids - 1 knee wedge - 1 positioning roll - 1 back-of-the head support - 2 large triangles - 2 small triangles - 1 wedge - 1 lordosis cushion - 2 flat wedges - 2 knee cushions - 2 different-sized multipurpose	 1 knee wedge  1 positioning roll




No.	Name of component	Description	Photograph
II.	Accessories:		
		cushions, 3 pieces each - Gap cushion set small - Gap cushion set large	  1 back-of-the-head support 2 large triangles   2 small triangles 1 wedge   1 lordosis cushion 2 flat wedges   2 knee cushions 2 different-sized multipurpose cushions, 3 pieces each   Gap cushion set small Gap cushion set large
26.	Paper Roller Holder Protective cover for table	Paper Roller Holder Protective cover for table	  1 paper roll holder Protective table cover


No.	Name of component	Description	Photograph
II.	Accessories:		
27.	Head Matrix Coil accessories for Head/Neck 16 Coil	Head Matrix Coil accessories - Relax cushion - Head/Neck cushion - Head stabilization wedge - Clamp cushions - Mirror carrier with mirror - filling cushion for feet first" scanning position.	 Relax cushion  Head/Neck cushion  Head stabilization wedges  Clamp cushions  Mirror carrier with mirror Head/Neck 16  Filling cushion for Feet First orientation
28.	Head Matrix Coil accessories for Head/Neck 20 Coil	Head Matrix Coil accessories - Relax cushion - Head/Neck cushion - Head stabilization wedge - Clamp cushions - Mirror carrier with mirror - Protective cover for the Head/Neck coil sockets	 Relax cushion  Head/Neck cushion  Head stabilization wedges  Clamp cushions




No.	Name of component	Description	Photograph
II.	Accessories:		
			  Mirror carrier with mirror Head/Neck 20 Protective cover for the Head/Neck 20 coil sockets
29.	TxRx CP Head accessories	TxRx CP Head accessories: - head cushion - head wedges (thin size) - head-spine table pad - head wedges (thick size)	  Head cushion 2 head wedges (thin size)   Head-Spine table pad 2 head wedges (thick size)
30.	Spine Coil accessories Spine 24 Coil Spine 32 Coil	Spine Coil accessories: - Cushion support	 Cushion support
31.	Body 6 Coil accessories	Body 6 Coil accessories: - 2 straps - 4 belts - iPAT cushion	  2 straps 4 belts



No.	Name of component	Description	Photograph
II.	Accessories:		
			 iPAT cushion
32.	Body 18 Coil accessories	Body 18 Coil accessories: - 2 straps - 4 belts - Cardiac spacer cushion	 2 straps  4 belts  Cardiac spacer cushion
33.	Accessories for Body coil for expert examinations Body 30/60 Coil	Body 30/60 accessories: - soft pad 1 for Body 30; - soft pad 2 for Body 30; - soft pad 3 for Body 30; - soft pad 4 for Body 30; - 2 belts; - 4 belts; - belt clamps; - frame; - adapter for Spine coil.	 Body 30 table pad 1  Body 30 table pad 2  Body 30 table pad 3  Body 30 table pad 4







No.	Name of component	Description	Photograph
II.	Accessories:		
			 2 straps  4 belts  4 strap clips  Frame for use as posterior coil  Adapter for spine coil sockets
34.	Accessories for Breast Coil for expert examinations Breast 18 coil	Accessories for Breast 18 coil: - coil pad; - body cushion 60; - body cushion 70; - head fixator; - headrest cushion; - arm cushion; - breast fixator cushion; - soft pad Breast 18.	 Coil pad  Body pad 60  Body pad 70  Head support  Head support cushion  Arm rest pad









No.	Name of component	Description	Photograph
II.	Accessories:		
			  Cushions for the breast fixations Table pad Breast 18
35.	Accessories for Biopsy coil 4-ch Breast BI Coil	Accessories for Biopsy coil: - Head rest - Abdominal pad - Foot support roller	   Head rest Abdominal pad Foot support roller
36.	Accessories for biopsy procedures for Modular Breast Coils for visualization and biopsy 2-/4-/8-ch Sentinelle Breast Coil 2-/10-/16-ch Sentinelle Breast Coil	Accessories for biopsy for 2-/4-/8 Sentinelle Breast Coil and 2-/10-/16 Sentinelle Breast Coil: - Fixed marker block/Marker block - Contra lateral support (for medial access) - Immobilization plate - Grid plate	   


No.	Name of component	Description	Photograph
II. Accessories:			
37.	Accessories for Modular Breast Coils for visualization and biopsy 2-/4-/8-ch Sentinelle Breast Coil 2-/10-/16-ch Sentinelle Breast Coil	Accessories for patient comfort for 2-/4-/8 Sentinelle Breast Coil and 2-/10-/16 Sentinelle Breast Coil: 10) Slider Breast Coil 11) Breast Cushion Set Body wedge pad; Shoulder wedges; Foot rest pad; Front pad; Hip wedge pad; Body tabletop pad; Insulation pad 12) Head rest 13) Head rest pad 14) Medial plug 15) Compression frame left 16) Immobilization plate 17) Phantom Holder 18) Compression frame right 19) Contralateral support 20) Patient support 21) Shoulder support; 22) Plastic Cover for 4 Ch coil element	
38.	Shallow Riser for Sentinelle Breast Coil	This option contains the base plate for 2-/4-/8-channel or 2-/10-/16-channel Sentinelle Breast Coil that needs to be positioned on the table.	
39.	Accessories for Breast Coil for visualization and biopsy Breast BI 7 Coil	Breast BI 7 accessories - Head rest - Abdominal pad - Foot support roller - Support cushion - Support cushion, narrow - Body pad 70 (hip wedge) - Armrest pad - Head support cushion - Fixation plates (lateral/medial)	  





No.	Name of component	Description	Photograph
II.	Accessories:		
		- Frame carrier lateral - Frame carrier medial - Cranio-caudal fixation unit (single element) - Pad for cranio-caudal fixation unit	 Support cushion  Support cushion, narrow  Body pad 70 (hip wedge)  Armrest pad  Head support cushion  Fixation plates (lateral/medial)  Frame carrier lateral  Frame carrier medial  Cranio-caudal fixation unit (single element)  Pad for cranio-caudal fixation unit
40.	Accessories for biopsy for Breast Coil for visualization and biopsy Breast BI 7 Coil	These accessories are intended for a biopsy for the coil for a biopsy and visualization of mammary glands. This set of accessories is intended only for training and not intended for patients Breast BI 7 accessories: - elevation plate (for medial access) - lateral grid - medial grid - reusable grid needle blocks (12G and 14G) - disposable grid needle blocks (12G, 14G, 16G and 18G)	 Elevation plate (for medial access)  Lateral grid  Medial grid  Reusable grid needle blocks (12G and 14G)  Disposable grid needle blocks (12G, 14G, 16G, and 18G)  Grid marker block




No.	Name of component	Description	Photograph
II.	Accessories:		
		- grid marker block - vertical post/pillar unit - needle guide for post/pillar unit - post/pillar marker - horizontal lateral slot plate - horizontal medial slot plate	 Vertical Post&Pillar unit  Needle guide for Post&Pillar unit  Post&Pillar marker  Horizontal lateral slot plate  Horizontal medial slot plate
41.	Accessories for Shoulder Coils Shoulder Small 16 Coil Shoulder Large 16 Coil	Shoulder Small/Large 16 accessories: - Base plate pad - Shoulder cup cushion	 Base plate pad  Shoulder cup cushion
42.	Accessories for Shoulder Shape 16 Coil	Shoulder Shape 16 Coil accessories: L-shaped holder	
43.	Accessories for Hand/Wrist 16 Coil	Hand/Wrist 16 accessories: - Wrist pad - Elbow pad - Palm pad - Wedge cushion - Wrist support cushion	 Wrist pad  Elbow pad  Palm pad  Wedge cushion






No.	Name of component	Description	Photograph
II.	Accessories:		
			 Wrist support cushion
44.	Accessories for TxRx CP Extremity Coil	TxRx CP Extremity coil accessories: - knee cushion; - leg support; - large foot ankle cushion; - small foot ankle cushion; - foot cushion; - air cushion; - foot support.	 Knee joint cushion  Leg support  Large ankle cushion  Small ankle cushion  Foot cushion  Air cushion  Foot support



No.	Name of component	Description	Photograph
II.	Accessories:		
45.	Accessories for TxRx Knee 15 coil TxRx Knee 15 Flare coil	Knee coils accessories: - Coil pad, 0.25 - Coil pad, 0.50 - Ankle pad - Knee rest pad - Foot pad	  Coil pad, 0.25 Coil pad, 0.50   Ankle pad Knee rest pad  Foot pad
46.	TxRx CP Extremity coil accessories	Knee coils accessories: - Coil pad, 0.25 - Coil pad, 0.50 - Ankle pad - Knee rest pad - Foot rest pad	  Coil pad, 0.25 Coil pad, 0.50   Ankle pad Knee rest pad




No.	Name of component	Description	Photograph
II.	Accessories:		
			 Foot rest pad
47.	TxRx Knee 18 coil accessories	Knee coils accessories: 5 – Knee cushion (knee inside coil) 3 – Leg cushion (the other leg) 4 – Feet cushion 2 – Cushion 4	
48.	Foot/Ankle 16 coil accessories	Foot/Ankle Coils accessories: - Knee cushion - Heel cushion - Foot cushion	 Knee cushion  Heel cushion  Foot cushion
49.	Coil accessories for Flex Coils Flex Large 4 Flex Small 4 UltraFlex Large 16 UltraFlex Small 16	Flex Coil accessories: - Hand wedge - 2 soft straps - 4 belts - iPat cushion, large - iPat cushion, small - 3 cable holders - Positioning pads for UltraFlex Coils - Soft wedge pad for UltraFlex coils - Foot/ankle positioning pads for UltraFlex coils	 Hand wedge  2 soft straps


No.	Name of component	Description	Photograph
II.	Accessories:		
			  4 belts iPat cushion, large   iPat cushion, small 3 cable holders  Positioning pads for UltraFlex coils  Soft wedge pad for UltraFlex coils  Foot/ankle positioning pads for UltraFlex coils

No.	Name of component	Description	Photograph
II.	Accessories:		
50.	Loop Coils accessories	Loop coil accessories: - Hand wedge - 2 soft straps - 4 belts - iPat cushion, large - iPat cushion, small - 3 cable holders	 Hand wedge  2 soft straps  4 belts  iPat cushion, large  iPat cushion, small  3 cable holders
51.	Peripheral Angio 36 Coil accessories	PA Matrix Coil accessories: - Leg support	 Leg support
52.	Endorectal Coil accessories	Endorectal Coil accessories - Endo cushion; - Endorectal Adapter 1.5T	 Endo cushion  Endorectal Adapter 1.5T


No.	Name of component	Description	Photograph
II. Accessories:			
53.	Prostate 2 coil accessories	Accessories for Prostate 2: - patient support pad; - back pad; - foot rest pad; - Body 18 cushion; - plug cap; - pressure testing cap; - cable extension; - base plate; - pressure testing unit	 Patient support pad  Back pad  Foot rest pad  Body 18 cushion  Plug cap  Pressure testing cap  Base plate  Pressure testing unit  Cable extension
54.	Accessories for Pediatric Coils	Accessories for pediatric coil: - body pad - head pad; - neck cushion; - head-spine pad; - head wedges; - Velcro straps;	 Body pad (for infant cradle only)  Head pad

No.	Name of component	Description	Photograph
II.	Accessories:		
			  Neck cushion Head-spine pad (for the use without infant cradle)   Head alignment wedges Velcro® straps (for infant cradle only)
55.	Vacuum cushions	Vacuum cushions for positioning and stabilization of the patient during the examination. The comfort kit may significantly reduce patient exam times and optimize image quality by minimizing the occurrence of motion artifacts caused by involuntary patient motion. The kit comprises a vacuum pump and 3 anatomically shaped cushions of different sizes (spine, head, and small joints).	  Head cushion Spine cushion  Small cushion
56.	Coil Storage Cart	Specially designed non-ferromagnetic cart for easy storage of the most commonly used coils and accessories.	

No.	Name of component	Description	Photograph
II. Accessories:			
57.	Optional Dot Engine Workflows	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.	
58.	Optional Application Packages	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.	 
59.	Optional Postprocessing Packages	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.	



No.	Name of component	Description	Photograph
II.	Accessories:		
60.	Additional optional single Applications	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.	 
61.	Argus Flow, Argus Dynamic Signal applications	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.	 
62.	syngo Expert-i, syngo Expert-i MRWP applications	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.	

No.	Name of component	Description	Photograph
II.	Accessories:		
63.	syngo System Security Enhanced, syngo System Security Enhanced applications	Software, optional application package. Software option for general regulatory security rules, providing enhanced security features including user management and audit trail functionality. This package supports customers in their achieving compliance with the HIPAA "Privacy" rule.	
64.	Video Distributor for monitor #T+D	Video distributor for connection of up to two additional color LCD monitors as well as an In-Room Acq WP syngo (in the examination room parallel to the syngo Acquisition Workplace operator monitor).	
65.	UPS system UPS Battery module	UPS system Liebert GXT4 3000RT230E for safeguarding computers. Including Power Cable of 9 m for connecting the UPS. UPS battery module Liebert GXT4 72VBATTE for safeguarding computers.	UPS 1ph 3kVA GXT4 + power cable  Battery module UPS GXT4 

No.	Name of component	Description	Photograph
II. Accessories:			
66.	Distribution board	Distribution board	
67.	RF-Cabin	Steel-full-room RF cabin, in a module method of building. Standard outside dimensions: L x W x H =6,590 x 4,090 x 2,750 mm Screen damping value: 90 dB (Co-Siting 100 dB) at 15 until 65 MHz	
68.	Filter plate for RF-Cabin	Filterplate in low leakage circuitry for shielded rooms (MRI) 6 separate lines, 250 V, 50/60 Hz, 250Vdc, 32A	
69.	Adapter, RF Cabin	2 adapter plates for adjusting the larger filter plate opening of an existing RF-room enclosure.	Filter plate with cabling: 



No.	Name of component	Description	Photograph
II. Accessories:			
70.	RF-Traverse Beams	5 traverse beams for the RF cabin.	Long Beam 6M, To be used with Lifting device  Lifting device without beams: 
71.	RF-Interior decor	Standard interior décor with Halogen/LED lighting for the RF cabin. Standard outside dimensions: L x W x H = 7,090 x 4,090 x 2,750 mm	
72.	RF Sound Proofing	Basic sound protection system. Standard outside dimensions of cabin: L x W x H = 7,090 x 4,090 x 2,750 mm.	
73.	Filter Cabinet, big Filter Cupboard	Cabinet for filter plate inside of the RF cabin. Wall-mounted cabinet custom-made. Outside dimensions appr. 1.300 x 420 x 2.700 mm (L x W x H) Cupboard for filter plate inside of the RF cabin.	



No.	Name of component	Description	Photograph
II.	Accessories:		
		Wall-mounted cabinet custom-made. Outside dimensions approx. 1,300 x 420 x 520 mm (L x W x H)	
74.	Threshold-free door	Soundproof threshold-free door for the RF-Cabins of MAGNETOM 1.5T/3T that can be opened and closed without effort by pneumatic mechanism. Clear opening: W x H 1,200 x 2,050 mm	
75.	Emergency hatch	With decorative plate on both sides, white, to be opened from the outside, approx. 600 x 600 mm.	
76.	Observation/ Day-light window	Observation / daylight window that can be opened, installation sideways. Dimensions: 1,200 x 800 mm.	
77.	Quench pipe	Quench pipe	

No.	Name of component	Description	Photograph
II. Accessories:			
78.	Eco Chiller 60kW	Eco Chiller 60kW The KKT ECO 133 -L chiller is a dedicated 20°C cooling system for all Tim 4G systems which automatically adapts to the different cooling requirements (e.g. system in operation, standby, ...) to reduce the energy consumption for cooling.	
79.	MR Workplace Table	The table design matches the MED-wide uniform design with silver-finished rim, use of friendly colors matching the Siemens color pattern for MAGNETOM and SOMATOM. - Width 120 cm - Depth 80 cm - Height 72 cm	 (1) MR Workplace Table (2) MR Workplace Container, 50 cm (optional)
80.	MR Workplace Container	The table design matches the MED-wide uniform design with silver-finished rim, use of friendly colors matching the Siemens color pattern for MAGNETOM and SOMATOM. Table height 72 cm, matching the syngo Acquisition Workplace and syngo MR Workplace console table, for installation in the operator room either directly to the left or right of the syngo Acquisition Workplace or syngo MR Workplace console table or separately. - Width 50 cm - Depth 80 cm	 (1) MR Workplace Table (2) MR Workplace Container, 50 cm (optional)

No.	Name of component	Description	Photograph
II.	Accessories:		
		- Height 72 cm Alternatively this casing is also suited for the Recon image processor (except for the MR systems with the new Tim generation: there the Recon image processor is always placed inside the electronics cabinet).	
81.	Glider Positioning Kit	Positioning System for patient transport between bed or wheelchair and MR-System. These sliding boards consist of a special material with very low coefficients of friction. They can be folded, take up little space and can be transported easily. They are firm but flexible and adapt to different shapes of the body. - Velcro safety belts fast and safe with Velcro closure. - Ergonomically shaped handles for safe and comfortable grasping - Radiolucent - Can be cleaned and disinfected with customary agents - Dimensions: 64 x 35 cm	

No.	Name of component	Description	Photograph
II.	Accessories:		
82.	Instrument table, non-magnetic	Mobile 4-leg frame construction made of 18/10 CrNi steel. Optimal - referring to the handle bar - working height of 800 mm. 63.5 cm x 43.5 cm tubular frame forms a safety rail; no sharp edges. Two shelves: surface ground to reduce glare, resistant to disinfectants, removable. Standard sizes, dimensioned to conform with the DIN standard for sterilized goods. 4 plastic swivel castors, two of them with brakes, wear resistant, non-discoloring when in contact with other material, electrically conductive, Ø 75 mm, 4 wall spacer rings.	
83.	Metal detector	Handheld detector for ferrous and nonferrous metals, non-contact examination of patients, battery or accu operated.	
84.	Fire Extinguisher non-magnetic	Non-magnetic fire extinguisher filled with carbon dioxide according to fire class B. With carbon dioxide filling, synthetic nozzle and wall holder for use in the MR suite.	

No.	Name of component	Description	Photograph
II.	Accessories:		
85.	Seismic Kit	Earthquake Kit for ECO Chiller 60 kW.	
86.	Service Syphon Kit	Service Syphon Kit	
	HE-SYPHON DEWAR END	HE-SYPHON DEWAR END	

No.	Name of component	Description	Photograph
II. Accessories:			
88.	SYPHON Extention	SYPHON EXT.LEG	
89.	Floor plate for table	Base Plate for the fixed table	
90.	Adhesive Plastic-Mastic complete kit	Adhesive Plastic-Mastic complete kit consisting of adhesive component A, adhesive component B, gloves and spatula.	
91.	Combi Dockable Table	Combi Dockable Table	Combi Dockable Table  (1) Connections (2) Tabletop with Transfer Board supplied by MAQUET (3) Poles MAQUET interface box 


No.	Name of component	Description	Photograph
II.	Accessories:		
			TRUMPF interface box  (1) Transfer board Maquet  Transfer Board TRUMPF  (2) Handle 

5.3 Technical requirements

Technical Requirements based on the key parameters and values as presented in the chapter below.

MagnaCoil™ Headset System



The MagnaCoil Headset System is an MR Conditional headphone. You can use it as an alternative headphone (for example, in head coils). The headset reduces the noise and can be used for communicating with the patient.

- (1) MagnaCoil Headset System incl. table plug and MagnaCoil Coupler for MagnaPlugs
- (2) MagnaPlugs (two sizes available: standard and mini)
- (3) MagnaCoil Coupler for MagnaPlugs

Intercom

The intercom allows personnel and patients to communicate during the examination. In addition, some important operations like stopping the patient table can be managed from the intercom. Optionally, music or automatic voice outputs can be played in the examination room via the loudspeaker or the headphones. The operating unit of the intercom is located at the syngo Acquisition Workplace.


- (1) Stop button
- (2) Speak button
- (3) Reset Patient Alert button
- (4) Listen button
- (5) Physio Signal button
- (6) Audio connection
- (7) Reset Stop button

Dot Display

The Dot (Day Optimized Throughput) display provides status information as well as several tools.

The display is located above the magnet bore opening at the front of the magnet cover. Optionally, at the back of the magnet, a second display can be installed for the additional control unit.


The display shows several pieces of information:

- Information concerning patient positioning (e.g. patient orientation, current table position, auto-positioning)
- Information about the connected coils
- Guidance for attaching ECG electrodes, docking the table and emergency evacuation
- General patient information submitted from the software
- Information concerning troubleshooting

Control Units

The control units are located at the right and the left side of the patient table at the front of the magnet cover. Optionally, additional control units are located at the back of the magnet.

Each control unit consists of a jogwheel and several additional buttons.

(1) Control unit for positioning the table

(2) Control unit for menu navigation and for adjusting settings for patient comfort

Depending on the situation, the small jogwheel can have different functions:

- Rotate: e.g. for navigating through the display menu or for changing settings of, for example, ventilation
- Push: e.g. for resetting the patient alert or for starting a measurement

5.3.1 Technical Requirements based on the key parameters and values as presented in the chapter below

The data and protocols included in this manual provide the operator with the means for evaluating the MR compatibility of non-Siemens MR products with the MR system described in the System Owner Manual. The MR compatibility of non-Siemens MR products addresses products of both manufacturers, the company producing the non-Siemens MR product and the company producing the MR system. The operator is ensured the safe operation of the non-Siemens MR product together with the MR system only when both manufacturers provide a test certificate addressing the "testing the effect of non-Siemens MR products on MR systems" as well as "testing the effect of MR systems on non-Siemens MR products". In all other cases, the operator has to ensure that use of the non-Siemens MR product does not conflict with the functions of the MR system and vice versa.

According to IEC 60601-2-33, the manufacturer of the MR system is obligated to provide a data sheet covering technical information of the MR system to enable items to be assessed for MR compatibility.

Own liability and risk

Based on the information provided, the operator evaluates the non-Siemens MR product at his own risk within the MR environment. Siemens shall not be held liable for any hazards resulting from this evaluation.

MR compatibility test

Manufacturers of a non-Siemens MR product, who would like to obtain a test certificate ("testing the effect of non-Siemens MR products on MR systems") for their device from Siemens, should get in touch with their Siemens sales engineer.

Parameters

All system-specific parameters are included in the "Technical Data" register of the System Owner Manual.

MR compatibility protocols

According to IEC 60601-2-33, protocols need to be proposed for testing the functionality of non-Siemens MR products. The protocols listed below are routinely used in the system.

The tests cannot be used to evaluate the effects of the non-Siemens MR product on the image quality of the MR system.

None of the service sequences listed below can be used to confirm MR compatibility or to define MR Conditional parameters of devices in the MR environment.

Program Card Sub-directory	Service Sequence	Parameters	Description & Comments
Default/ Default/ ServiceSequences/ Default	se_multiecho	Use default parameters or change TR on the Routine parameter card to minimum or according to SAR limits to get max. RF power load. On the Routine parameter card change FoV and slice thickness to minimum to get max. gradients.	Multiple RF pulses and gradient pulses.
	rf_field	Use default parameters or change TR on the Routine parameter card to minimum or according to SAR limits to get max. RF power load. On the Routine parameter card change FoV and slice thickness to minimum to get max. gradients.	Long (4 ms) powerful RF pulse.
	tuncal	Use default parameters or change Flip angle on the Contrast parameter card to maximum according to SAR limits to get max. RF power load. On the Resolution parameter card change Phase resolution to maximum to increase scan time.	RF power only. Repetitive RF pulses. The system-specific frequency band is measured in steps.
	spike	Use default parameters. In case "gradient stimulation limit exceeded" occurs change orientation to coronal or sagittal or run sequence under Service_Patient.	Gradients only. Repetitive multiple strong rapid gradient pulses.

Magnetic fringe field and control area

This table shows the effects of the magnetic fringe field on devices located in the vicinity of the magnet and the safety distances required. Observe the minimum distances to be maintained from the center of the x, y, and z axes of the magnet.

1 mT = 10 Gauss

Magnetic flux density	Minimum distances (x = y = radial, z = axial)	Examples: Devices affected
3 mT	x = 1.84 m z = 2.79 m	Small motors, watches, cameras, credit cards, magnetic media
1 mT	x = 2.20 m z = 3.46 m	Oscilloscopes, computers, disk drives, shielded color monitors
0.5 mT	x = 2.50 m z = 4.00 m	B/W monitors, magnetic media, cardiac pacemakers, insulin pumps
0.2 mT	x = 2.96 m z = 4.88 m	Siemens CT systems
0.1 mT	x = 3.1 m z = 5.5 m	Siemens linear accelerators
0.05 mT	x = 3.9 m z = 6.8 m	X-ray I.I., gamma cameras, third party linear accelerators

Static magnetic field B₀

The figures show lines of the same magnetic flux density in milliTesla.

The following graphics show the calculated magnetic field in air. Magnetic materials in the vicinity of the magnet (i.e. iron beams or room shielding) may influence the form of the stray field.

The plots represent three orthogonal planes through the isocenter to illustrate maximum spatial extent of iso-magnetic contours.

Each plot contains the iso-magnetic contours with values of 0.5 mT, 1 mT, 3mT, 5 mT, 10 mT, 20 mT, 40 mT, and 200 mT as well as a distance scale and a superimposed outline of the MR system.


The 0.5 mT line marks the exclusion zone of the static magnetic field (pace maker limit).


The magnetic field is rotationally symmetric about the z-axis and mirror symmetric about the horizontal axis.

The multiple plots are shown to provide information about the accessible space in different orientations. This note also applies to the other magnetic field related plots.

View in the direction of the magnet axis


Side view of the magnet

(1) Z Axis (m)


(2) Vertical Axis (m)

Top view of the magnet


Iso-magnetic contours at positions accessible to and relevant for MR workers

In the following, we show plots representing the 0.5 T, 1 T, and 1.5 T isomagnetic contours at positions accessible to and relevant for the MR worker as far as the static magnetic field in the isocenter exceeds any of these values. The plots all represent one quarter of the magnet as indicated below. The plot lines each run symmetric to the axis to the isocenter of the magnet. As a result, the graphs can be unfolded.

Top view


- (1) Z Axis (m)
- (2) X Axis (m)
- (3) Magnet

Side view

- (1) Z Axis (m)
- (2) Y Axis (m)
- (3) Patient table

Side view


- (1) Z Axis (m)
- (2) Vertical (m)
- (3) Patient table

⊗: At this location, the value of the magnetic field B_0 is greatest.

The iso-magnetic contour lines are accurate to a value better than 1%. There is very little influence from the environment to the scanner. The graphic shows the cover in nominal position. On a variety of scanners, the position can vary ± 5 mm in the axial direction. The maximum value therefore has a tolerance of $\pm 10\%$, as small geometric deviations can cause a significant change in the value.

Top view

(1) Z Axis (m)

(2) Horizontal (m)

⊗: At this location, the value of the magnetic field B_0 is greatest.


Spatial gradient of the static magnetic field B_0

The rise of the magnetic field as a function of the distance to the magnet is expressed by the spatial gradient of B_0 . The following figures show lines with the same gradient in T/m. The magnetic attraction force on a magnetically saturated ferromagnetic object is proportional to this quantity.

Please note: Sometimes this quantity is expressed in G/cm (1 T/m = 100 G/cm).

Side view


- (1) Z Axis (m)
- (2) Vertical (m)
- (3) Patient table

⊗: At this location, the force on a magnetically saturated ferromagnetic object is greatest.

The iso-magnetic contour lines are accurate to a value better than 1%. There is very little influence from the environment to the scanner. The graphic shows the cover in nominal position. On a variety of scanners, the position can vary ± 5 mm in the axial direction. The maximum value therefore has a tolerance of $\pm 10\%$, as small geometric deviations can cause a significant change in the value.

Top view

(1) Z Axis (m)


(2) Horizontal (m)

⊗: At this location, the force on a magnetically saturated ferromagnetic object is greatest.

Product of the static magnetic field B_0 and the spatial gradient of B_0

The magnetic attraction force on a diamagnetic/paramagnetic object or on a ferromagnetic material below its magnetic saturation is proportional to this quantity.

Side view


- (1) Z Axis (m)
- (2) Vertical (m)
- (3) Patient table

⊗: At this location, the force on a diamagnetic/paramagnetic object or on a ferromagnetic material below its magnetic saturation is greatest.

The iso-magnetic contour lines are accurate to a value better than 1%. There is very little influence from the environment to the scanner. The graphic shows the cover in nominal position. On a variety of scanners, the position can vary ± 5 mm in the axial direction. The maximum value therefore has a tolerance of $\pm 10\%$, as small geometric deviations can cause a significant change in the value.

Top view

(1) Z Axis (m)

(2) Horizontal (m)

⊗: At this location, the force on a diamagnetic/paramagnetic object or on a ferromagnetic material below its magnetic saturation is greatest.

Gradient stray field distribution relevant for assessing exposure to MR workers

The diagrams show the magnetic stray field distribution of the gradient system along the patient axis of the MAGNETOM Aera with the XJ or XQ gradient system respectively according to the requirements of the IEC 60601-2-33 standard.

The magnetic fields generated by each of the three orthogonal gradient axes are calculated on a coordinate grid on a virtual cylinder surface. The cylinder encompasses the patient axis and starts at the magnet isocenter.


The grid is defined by:

- a point pattern on the surface of a virtual cylinder with a diameter of 70 cm (i.e. max. bore liner);
- points on circles on the cylinder surface perpendicular to the cylinder axis with an azimuthal step width of 10 degrees (i.e. 36 points at the circumference) and an axial step width of 4 cm.

The field values are calculated at max. currents that can be handled by the gradient amplifiers. At each space point the field vectors generated by each gradient are calculated, superimposed and then their magnitudes are derived.


The maximum magnitude field value that can be found on each circle is selected and plotted along the z-axis. By dividing the maximum field values by the shortest rise time, the dB/dt values can be derived.

	Shortest rise time	Max dB/dt
XJ gradient	264 μ s	57 mT/264 μ s = 216 T/s
XQ gradient	225 μ s	70 mT/225 μ s = 311 T/s

XJ gradient system

- (1) Magnet isocenter
- (2) End of magnet

XQ gradient system


- (1) Magnet isocenter
- (2) End of magnet

The areas shaded in grey indicate the length of the magnet.

Spatial distribution of the gradient field at different positions inside the gradient coil


The diagrams show the spatial distribution of the maximum magnitude values of the vector sum of the field components generated by each of the three gradient units simultaneously at positions on virtual cylinders coaxial with the patient axis with diameters of 0.2 m, 0.4 m, and 0.6 m.

The virtual cylinders have the same length as the gradient coil. In the axis direction of the cylinder, the points have a separation of 0.02 cm.

XJ gradient system

- (1) Magnet isocenter
- (2) End of gradient coil


XQ gradient system

- (1) Magnet isocenter
- (2) End of gradient coil

The areas shaded in grey indicate the length of the gradient coil.

RF power distribution relevant for assessing exposure to MR workers

The diagram shows the RF power distribution along the patient axis of the MR system according to the requirements of the IEC 60601-2-33 standard.

The RF field is calculated along the center line starting in the magnet isocenter.

The area shaded in grey indicates the length of the system.

The ratio $B1(z)^2/B1(0)^2$ provides a worst case estimation of the SAR contribution to a person who is positioned at a distance z from the isocenter. The SAR contribution is relative to the SAR that is applied to a person in the center of the patient bore.

For example a person standing in front of the system aperture absorbs a maximum of 0.2% of the RF power which is applied to a patient scanned in the center of the bore.

5.3.2 Guidance and manufacturer's declaration EMC

Preamble

The following document applies to all products provided by the Magnetic Resonance (MR) product group of Siemens Healthineers.

Deviations and additions to this document are provided in accompanying product-specific documents.

Basic safety and ESSENTIAL PERFORMANCE

As a result of the risk management process for the MR systems and according to IEC 60601-1: 2012 and IEC 60601-2-33:2015, no ESSENTIAL PERFORMANCE was identified.

MR systems provide basic safety and ESSENTIAL PERFORMANCE in the electromagnetic environments of the intended use throughout their expected service life. Please refer to the "Maintenance Plan" for all instructions necessary to maintain basic safety and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the expected service life.

Environment, RF shielding and RF interferences

Environment

For the MR system, two environmental categories according to IEC 60601-1-2:2014 are defined:

- Professional healthcare facility environment

This environment covers the control room with the MR console/host as well as the equipment room where all hardware equipment outside of the RF-shielded examination room is installed.

- Special environment: RF-shielded room

This environment covers the examination room where the magnet and all hardware equipment inside of the RF-shielded room are installed.

The MR system in an RF-shielded room is suitable for use in professional healthcare facility environment, such as clinics and hospitals.

RF shielding

An RF-shielded room has to meet the following specifications:

7 T, including multi-nuclei spectroscopy option:

- frequency range: 13.5 MHz–298 MHz
- minimum RF shielding effectiveness: > 65 dB
- RF attenuation of the RF cabin itself (without any feedthroughs or filters): > 90 dB

1.5 T and 3 T, including multi-nuclei spectroscopy option:

- frequency range: 16.7 MHz–123.2 MHz
- minimum RF shielding effectiveness: > 65 dB
- RF attenuation of the RF cabin itself (without any feedthroughs or filters): > 90 dB

It is essential that the actual shielding effectiveness and filter attenuation of the shielded location be verified to assure that they meet the minimum specification. Verification of the shielding requirements is essential prior to handover of the system. The integrity of the RF shielding is the responsibility of the customer.

Recommended test methods for measurement of RF shielding effectiveness and RF filter attenuation are given in the "Planning Guide System".

RF operation of MR system

MR systems apply RF energy for diagnosis. Operation involving RF energy is performed in the RF-shielded room. Frequency bands for operation of MR systems are as follows.

63.6 MHz	1.5 T
16.7 MHz–123.2 MHz	3 T, including multi-nuclei spectroscopy option
13.5 MHz–298 MHz	7 T, including multi-nuclei spectroscopy option

The MR system applies IQ-modulated (without encoding any data) repeated RF pulses. This results in a total bandwidth of the modulated signal of less than ± 500 kHz around the carrier frequency.

Compatible equipment

In principle, MR systems may be used in combination with compatible equipment only. For details on compatibility, see: Operator Manual MR System.

The use of the MR system adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, the MR system and the other equipment should be observed in order to verify that they are operating normally.

Please note, HF surgical equipment is not compatible with MR systems.

Accessories and cabling

All cables, transducers, and accessories are part of the MR system and cannot be removed or replaced by the user. Unauthorized replacement of accessories and cabling may affect the EMC behavior of the MR system and cause deviation from the EMC requirements of IEC 60601-1-2:2014, and furthermore may influence the functionality of the system.

The use of accessories, transducers, and cables other than those specified, with the exception of transducers and cables sold by the manufacturer of the equipment or system as replacement parts for internal components, may result in increased electromagnetic emission or decreased electromagnetic immunity of the equipment or system and result in improper operation.

RF interference

MR systems apply and receive RF energy for diagnosis. Thus, MR systems may interfere with other equipment located nearby.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the MR system, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Failure to use the MR system in the specified type of shield location could result in degradation of the performance of this equipment, interference with other equipment, or interference with radio services.

The use of RF emitting equipment, which has not been released or certified by Siemens, inside the shielded location of the MR system is not allowed unless the equipment is switched to an RF silent mode before being placed in the shielded location. In case of doubt, please check the operator manual or the guidelines of the respective third-party equipment.

Electromagnetic emission

According to IEC 60601-2-33:2015, the MR system outside the special environment is in compliance with emission requirements of clause 7 of

IEC 60601-1-2:2014. Inside the special environment, clause 7 of IEC 60601-1-2:2014 does not apply.

Emission tests	Compliance
Radiated RF emission (150 kHz–1 GHz)	EN55011/CISPR 11 Class A/Group 2
Conducted RF emission (150 kHz–30 MHz)	EN55011/CISPR 11 Class A/Group 2

The emissions characteristics of this equipment make it suitable for use in a professional healthcare facility environment, such as clinics and hospitals (CISPR 11 class A). If it is used in a residential environment, this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

Electromagnetic immunity

General requirements

According to IEC 60601-2-33:2015, the MR system is intended for use in the electromagnetic environment specified as follows:

- The MR system outside the special environment complies with clause 8 of IEC 60601-1-2:2014.
- All operator and patient accessible parts in the special environment are compliant with the ESD requirements in clause 8.1 of IEC 60601-1-2:2014.

The customer or the user of the MR system should assure that it is used in such an environment. Normally, the MR system is used in a professional healthcare facility environment.

Immunity tests	Test level and compliance level
Electrostatic discharge, air discharge	EN/IEC 61000-4-2 ± 2.0 kV, ± 4.0 kV, ± 8.0 kV, ± 15.0 kV
Electrostatic discharge, contact discharge	EN/IEC 61000-4-2 ± 8.0 kV
Radiated RF electromagnetic fields	EN/IEC 61000-4-3 80 MHz–2.7 GHz, 3 V/m, 1 kHz, 80 % AM

Immunity tests	Test level and compliance level
Proximity fields from RF wireless communications equipment	EN/IEC 61000-4-3 385 MHz: 27 V/m 450 MHz: 28 V/m 710 MHz/745 MHz/780 MHz: 9 V/m 810 MHz/870 MHz/930 MHz: 28 V/m 1720 MHz/1845 MHz/1970 MHz: 28 V/m 2450 MHz: 28 V/m 5240 MHz/5500 MHz/5785 MHz: 9 V/m
Burst Power supply lines	EN/IEC 61000-4-4 ± 2.0 kV, 100 kHz repetition frequency
Burst Signal input/output ports	EN/IEC 61000-4-4 ± 1.0 kV, 100 kHz repetition frequency
Surge Power supply lines (line to ground)	EN/IEC 61000-4-5 ± 0.5 kV, ± 1.0 kV, ± 2.0 kV at 0°, 90°, 180° and 270°
Surge Power supply lines (line to line)	EN/IEC 61000-4-5 ± 0.5 kV, ± 1.0 kV at 0°, 90°, 180° and 270°
Conducted disturbances, induced by RF fields	EN/IEC 61000-4-6 1 kHz, 80 % AM 3 V, 150 kHz–80 MHz 6 V, in ISM and amateur radio bands between 150 kHz and 80 MHz
Magnetic fields	EN/IEC 61000-4-8 30 A/m, 50 Hz
Voltage interruptions on power supply input lines	EN/IEC 61000-4-11 Voltage interruptions: 0% UT 250/300 cycles

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

In situ test

If the exemption specified in IEC 60601-1-2, section 8.6 was applied for a permanently installed large MR system, the MR system may be tested in situ with the general requirements defined above except for the radiated RF immunity over the entire frequency range from 80 MHz to 6 GHz.

The MR system has been tested for radiated RF immunity only at selected frequencies, and use nearby of emitters at other frequencies could result in improper operation.

5.4 Key parameters & values provided by manufacturer

5.4.1 syngo MR E11E

Magnet System

Magnet parameters	
Operating field strength	1.5 Tesla
Magnet type	Superconductor
Field stability over time	<0.1 ppm/h
Weight (with cryogens)	3118 kg
Magnet length	137 cm
System length cover to cover	145 cm
Open Bore design ¹	70 cm

1. incl. shim coils, gradient coil, RF body coil

Homogeneity - MAGNETOM Aera XJ gradients:

Homogeneity (based on highly accurate 24 plane plot)		
	Guaranteed	Typical
10 cm DSV	0.02 ppm	0.01 ppm
20 cm DSV	0.075 ppm	0.06 ppm
30 cm DSV	0.3 ppm	0.2 ppm
40 cm DSV	1.4 ppm	1.1 ppm
50x50x45 cm ³ DEV	4.0 ppm	3.1 ppm

In compliance with the German "Qualifikationsvereinbarung".

Standard deviation Vrms (Volume root-mean square) measured with highly accurate 24 plane plot method (20 points per plane).

Standard active shim with 3 linear channels.

DSV – Diameter spherical volume (x, y and z direction).

DEV – Diameter elliptical volume.

Homogeneity - MAGNETOM Aera XQ gradients:

Homogeneity (based on highly accurate 24 plane plot)		
	Guaranteed	Typical
10 cm DSV	0.01 ppm	0.003 ppm
20 cm DSV	0.05 ppm	0.03 ppm
30 cm DSV	0.3 ppm	0.15 ppm
40 cm DSV	1.4 ppm	1.0 ppm
50x50x45 cm ³ DEV	4.0 ppm	3.0 ppm

In compliance with the German "Qualifikationsvereinbarung".

Standard deviation Vrms (Volume root-mean square) measured with highly accurate 24 plane plot method (20 points per plane).

Active shim with 3 linear and 5 non-linear channels (2nd order).

DSV – Diameter spherical volume (x, y and z direction).

DEV – Diameter elliptical volume.

Shimming	
Both: passive and active shimming. Passive shimming during installation.	
Standard active shim with 3 linear channels (1 st order). Additional 5 non linear channels (2 nd order).	
3D Shim	Patient-specific automated shim
	Time to shim = approx. 20 s

Shielding		
Active Shielding (AS)	5 th generation active shielding (AS) technology with counter coils	
Fringe field (axialxradial)	0.5 mT ¹	4.0 m×2.5 m
	0.1 mT	5.5 m×3.1 m
External Interference Shield (E.I.S.)	Patented shielding system integrated into the magnet	
	Continuous compensation and automatic suppression of external magnetic field interferences during measurement (caused by moving ferromagnetic objects or nearby power lines)	

1. pacemaker safety limit

Magnet cooling system	
Refill interval (typical) ¹	Not applicable
Boil-off rate (typical) ¹	0.0 l/year
Max. helium capacity	approx. 1280 liters
Minimum helium level	35%
Cryostat	Stainless steel

1. For typical clinical use, depending on sequences and operating time with running helium compressor. The system needs to be serviced at regular interval. Undisturbed magnet cooling for 24 hours and 7 days a week.

Patient handling

The patient table is available in two different configurations:

Tim table		
Max. patient weight for vertical and horizontal table movement		250 kg
Max. scan range		140 cm optional 205 cm ¹
Vertical table movement	Range	52–102 cm ² +13 mm ³
	Speed	6 cm/s
Horizontal table movement	Max. range	261 cm
	Max. speed	20 cm/s
	Horizontal accuracy for repositioning from one direction	±0.5 mm
Continuous table movement during scan capable		

Tim Dockable Table		
Max. patient weight for vertical and horizontal table movement		250 kg
Max. scan range		140 cm optional 205 cm ¹
Vertical table movement	Range	56–106 cm ² +13 mm ³
	Speed	6 cm/s
Horizontal table movement	Max. range	261 cm
	Max. speed	20 cm/s
	Horizontal accuracy for repositioning from one direction	±0.5 mm
Continuous table movement during scan capable		

Hearing protection data

The A-evaluated, effective sound pressure level was measured according to NEMA MS 4-2006 (National Electrical Manufacturers Association) using the maximum gradient acoustic noise (MGAN) method.

Patient noise - MAGNETOM Aera XJ gradients:

Patient noise
Patients require hearing protection with an SNR ¹ = 18 dB or more

1. SNR = Single Number Rating

Patient noise - MAGNETOM Aera XQ gradients:

Patient noise
Patients require hearing protection with an SNR ¹ = 22 dB or more

1. SNR = Single Number Rating

Noise for personnel in the examination room - MAGNETOM Aera XJ gradients:

Noise for personnel in the examination room	
Noise measured	83.1 dB(A)
Hearing protection	SNR = 1 dB or more

Noise for personnel in the examination room - MAGNETOM Aera XQ gradients:

Noise for personnel in the examination room	
Noise measured	86.1 dB(A)
Hearing protection	SNR = 4 dB or more

The standard Siemens headphones are intended for communication with the patient and can be used in combination with ear plugs.

Gradient system

MAGNETOM Aera XJ gradients:

Gradient amplifier	
• Water-cooled, highly compact, modular design • Ultra-fast solid-state technology with very low switching losses	
Max. output voltage ¹	2000 V
Max. output current ¹	625 A

1. Values for each of the 3 gradient axes

MAGNETOM Aera XQ gradients:

Gradient amplifier	
• Water-cooled, highly compact, modular design • Ultra-fast solid-state technology with very low switching losses	
Max. output voltage ¹	2250 V
Max. output current ¹	900 A

1. Values for each of the 3 gradient axes

Gradient System Performance- MAGNETOM Aera XJ gradients:

Performance	
Performance for each axis	
Max. amplitude	33 mT/m
Min. rise time	264 μ s
Max. slew rate	125 T/m/s
Vector gradient performance (vector addition of all 3 gradient axes)	
Max. eff. amplitude	57 mT/m
Max. eff. slew rate	216 T/m/s
Gradient duty cycle	100%

Gradient System Performance – MAGNETOM Aera XQ gradients:

Performance	
Performance for each axis	
Max. amplitude	45 mT/m
Min. rise time	225 μ s
Max. slew rate	200 T/m/s
Vector gradient performance (vector addition of all 3 gradient axes)	
Max. eff. amplitude	78 mT/m
Max. eff. slew rate	346 T/m/s
Gradient duty cycle	100%

Direct transmit technology		
Frequency stability (5 min)	$\pm 2 \times 10^{-10}$	
Frequency control	32 bits (0.015 Hz)	
Phase control	16 bits (0.006 degrees)	
Body coil	Integrated whole body no tune transmit/receive coil with 16 rungs	
	Optimized RF efficiency and signal-to-noise ratio (SNR)	
Transmitter path	Feedback loop for excellent RF stabilization	
	Transmit amplitude	16 bit control 25 ns resolution
	Gain stability (after first minute)	<0.05 dB (1 s) <0.2 dB (5 min)
Transmitter amplifier	Extremely compact, water-cooled solid state amplifier, integrated at the magnet as part of DirectRF™ technology	
	Transmit amplifier bandwidth	800 kHz
	Peak power	26.1 kW

RF receiver technology	
Number of coil elements	Up to 204
Number of independent receiver channels	24, 48, 64 ¹
Quadrature demodulation and filtering	Digital
Receiver bandwidth	500 Hz–1 MHz (for each channel)
Receiver signal resolution	32 bit
ADC sampling rate	80 MHz
Preamplifier noise figure	<0.5 dB
Dynamic range at coil connector (referred to 1 Hz resolution bandwidth)	164 dB instantaneous at receiver 169 dB with automatic gain control at local coil connector

1. Optional

syngo Acquisition Workplace		
Host computer	Processor	Intel Xeon ≥ E5-1620 Quad-Core
	Clock rate	3.6 GHz
	Main memory (RAM)	32 GB
	1 st hard disk (system SW)	≥300 GB SAS
	2 nd hard disk (data base)	≥300 GB SAS
	3 rd hard disk (images)	≥300 GB SAS
	CD-R writer	Approx. 4000 images 256 ² ; DICOM Standard, ISO 9660
	DVD-R writer	Approx. 25000 images 256 ² ; DICOM Standard, ISO 9660
	Media drives	CD/DVD drive
Color LCD monitor ¹	Screen size (diagonal)	19"
	Horizontal frequency	30–100 kHz
	Vertical frequency	50–75 Hz
	Screen matrix	1280×1024 pixels
Measurement and reconstruction system (system configuration with 24 receiver channels)	Processor	Intel ≥ E3-1225 Quad-Core
	Clock rate	3.2 GHz, or comparable
	Main memory (RAM)	32 GB
	Hard disk for raw data	≥300 GB
	Hard disk for system software	≥100 GB
	Parallel Scan and Recon	Simultaneous scan and reconstruction of up to 12 data sets
	Reconstruction speed	12987 recons per second (256 ² FFT, full FoV) 57971 recons per second (256 ² FFT, 25% recFoV)

1. A standard monitor without calibration is not suitable for diagnostic purposes. Please consider the initial acceptance testing for image display devices and the follow-up service for constancy testing on a regular base, as offered by Siemens service.

syngo Acquisition Workplace		
Measurement and reconstruction system (system configuration with 48 receiver channels)	Processor	Intel $\geq$ E5-2620 2x6-Core
	Clock rate	2x2.0 GHz, or comparable
	Main memory (RAM)	48 GB
	Hard disk for raw data	$\geq$ 300 GB
	Hard disk for system software	$\geq$ 100 GB
	Parallel Scan and Recon	Simultaneous scan and reconstruction of up to 12 data sets
	Reconstruction speed	22556 recons per second (256 ² FFT, full FoV) 88889 recons per second (256 ² FFT, 25% recFoV)

syngo Acquisition Workplace		
Measurement and reconstruction system (optional system configuration with 64 receiver channels)	Processor	Intel Xeon E5-2658 2x8-Core
	Clock rate	$\geq$ 2x2.1 GHz
	Main memory (RAM)	64 GB
	Hard disk for raw data	$\geq$ 400 GB
	Hard disk for system software	$\geq$ 100 GB
	Reconstruction speed	31128 recons per second (256 ² FFT, full FoV) 118519 recons per second (256 ² FFT, 25% recFoV)
	Parallel Scan and Recon	Simultaneous scan and reconstruction of up to 12 data sets
	GP GPU	1xTesla K10

syngo MR Workplace (optional)

Color LCD monitor and host computer as for syngo Acquisition Workplace

Euro Connector (System IEC 320)

Load rating	100–240 Vac
	50–60 Hz
	max. 10 A

Cooling system

Two different customer specific cooling alternatives (Separator or Eco Chiller) are available.

Cooling system options – MAGNETOM Aera XQ gradients:

Separator option (for connection to available cooling system)	Water consumption	70 l/min ¹
	Heat dissipation to water	45 kW
Eco Chiller option with automatic adaptation to the required cooling demands (e.g. different night/day mode) to decrease energy cost	GREEN Cooling Package²: • Reduces energy consumption by up to 50% • Automatic start if the surrounding temperature is 18 °C (64 °F) or less • If the temperature is less than –3 °C (27 °F) the chiller is switched off³	

1. Water temperature: 12 °C (54 °F)

2. Free Cooling Unit, optional

3. In case of clinical routine measurement conditions

Cooling system options – MAGNETOM Aera XQ gradients:

Separator option (for connection to available cooling system)	Water consumption	90 l/min ¹
	Heat dissipation to water	60 kW
Eco Chiller option with automatic adaptation to the required cooling demands (e.g. different night/day mode) to decrease energy cost	GREEN Cooling Package²: • Reduces energy consumption by up to 50% • Automatic start if the surrounding temperature is 18 °C (64 °F) or less • If the temperature is less than –10 °C (14 °F) the chiller is switched off³	

1. Water temperature: 12 °C (54 °F)

2. Free Cooling Unit, optional

3. In case of clinical routine measurement conditions

Line power supply

	Values	Tolerance
Voltage	380 V, 400 V, 420 V, 440 V, 460 V, 480 V	±10%
Frequency	50/60 Hz	±1 Hz
Connection value	63 kVA for 480 V/400 V, 62 kVA for 380 V	

For all products, line power has to be supplied via an on-site system contact or via another multipole shut-down mechanism. Room installation has to be in compliance with VDE 0100-710¹.

1. In all countries, compliance with local and national legal regulations is required. However, we strongly recommend compliance with the regulations described herein – to the extent permitted by relevant local and national laws – in order to ensure the safety of operating personnel, patients, and third parties.

Power rating plate - MAGNETOM Aera XJ gradients

380/400/420/440/460/480	V ~
3PHASE	50/60 Hz
SHORT – TIME	kVA: 75
LONG – TIME	kVA: 20

Power rating plate - MAGNETOM Aera XQ gradients

380/400/420/440/460/480	V ~
3PHASE	50/60 Hz
SHORT – TIME	kVA: 102
LONG – TIME	kVA: 20

Power consumption- MAGNETOM Aera XJ gradients

Off ¹	6.1 kW
System ready to measure ¹	8.7 kW
Scan ¹	18.9 kW
Remarks: • All values are typical values, applicable for 400 V/50 Hz. • Consumption for optional separator pump not included.	

1. The power consumption described herein is based on results that were achieved in a setting according to the COCIR methodology MRI - Measurement of the energy consumption (<http://www.cocir.org/site/index.php?id=46>). Since many variables impact power consumption (e.g. sequences used for scanning and sequence parameters, scan time), there can be no guarantee that each customer will achieve the same values.

Power consumption- MAGNETOM Aera XQ gradients

Off ¹	6.1 kW
System ready to measure ¹	9.3 kW
Scan ¹	22.5 kW
Remarks: • All values are typical values, applicable for 400 V/50 Hz. • Consumption for optional separator pump not included.	

1. The power consumption described herein is based on results that were achieved in a setting according to the COCIR methodology MRI - Measurement of the energy consumption (<http://www.cocir.org/site/index.php?id=46>). Since many variables impact power consumption (e.g. sequences used for scanning and sequence parameters, scan time), there can be no guarantee that each customer will achieve the same values.

Dimensions

Component	Width [cm]	Depth [cm]	Height [cm]	Weight [kg]	Heat dissipa- tion [kW]
Examination Room					
Magnet 1.5 Tesla AS (incl. Helium)	205	137	215	3118	
Magnet in operation, incl. gradient coil, body coil, Tim table, and covers	231	405 433 ¹	219	4798	
Tim table	76	249	52-102 ² +13 mm ³		
Required min. room height clearance			240 ⁴		
Min. transport dimen- sions	231	155	214		
Control Room					
syngo Acquisition Workplace (table+monitor)	120	80	117 (72+45)		
Host computer	22	46	47		
syngo MR Workplace (optional) (table+monitor)	120	80	117 (72+45)		
Equipment Room					
Electronics cabinet, incl. system control, RF system, gradient power system, image processor	160	65	198 ⁵	1500	≤5 ⁶
Cooling system	65	65	189	500	

Additional system characteristics

Magnet data	
Operating field strength	1.5 Tesla
Magnet type	Superconductor
Field stability over time	<0.1 ppm/h
Weight (with cryogenics)	3118 kg
Magnet length	137 cm
System length cover to cover	145 cm
Open Bore design ¹	70 cm
Refill interval (typical) ²	Not applicable
Boil-off rate (typical) ²	0.0 l/year
Max. helium capacity	approx. 1280 liters
Minimum helium level	35%
Cryostat	Stainless steel

RF data		
Peak power of transmitter amplifier	26.1 kW	
Transmitter bandwidth	800 kHz	
Receiver bandwidth	500 Hz–1 MHz (for each channel)	
RF transmit coils		
Body coil	Max. applied RF field B1+	27 μ T
	Max. specified B1+ rms	6.2 μ T
	Distance to isocenter, where RF transmit field is reduced by	
	3 dB	0.18 m
	10 dB	0.38 m
CP Extremity	Max. applied RF field B1+	47 μ T
	Max. specified B1+ rms	17 μ T
TxRx 15 Ch Knee	Max. applied RF field B1+	47 μ T
	Max. specified B1+ rms	13 μ T
TxRx Knee 15	Max. applied RF field B1+	47 μ T
	Max. specified B1+ rms	13 μ T
TxRx CP Head	Max. applied RF field B1+	47 μ T
	Max. specified B1+ rms	12.5 μ T

Tim table		
Max. patient weight for vertical and horizontal table movement		250 kg
Max. scan range		140 cm; optional 205 cm ¹
Vertical table movement	Range	52–102 cm ² +13 mm ³
	Speed	6 cm/s
Horizontal table movement	Max. range	261 cm
	Max. speed	20 cm/s
	Horizontal accuracy for repositioning from one direction	±0.5 mm
Continuous table movement during scan capable		

Tim Dockable Table		
Max. patient weight for vertical and horizontal table movement		250 kg
Max. scan range		140 cm; optional 205 cm ¹
Vertical table movement	Range	56–106 cm ² +13 mm ³
	Speed	6 cm/s
Horizontal table movement	Max. range	261 cm
	Max. speed	20 cm/s
	Horizontal accuracy for repositioning from one direction	±0.5 mm
Continuous table movement during scan capable		

Patient comfort	
Open bore design ¹	70 cm
In-bore lighting	can be set at 6 different levels
In-bore ventilation	can be set at 6 different levels
In-bore intercom	including loudspeaker, microphone and earphones

5.4.2 syngo MR E11C

Magnet System

Magnet parameters	
Operating field strength	1.5 Tesla
Magnet type	Superconductor
Field stability over time	<0.1 ppm/h
Weight (with cryogenics)	3118 kg
Magnet length	137 cm
System length cover to cover	145 cm
Open Bore design ¹	70 cm

Homogeneity - MAGNETOM Aera XJ gradients:

Homogeneity (based on highly accurate 24 plane plot)		
	Guaranteed	Typical
10 cm DSV	0.02 ppm	0.01 ppm
20 cm DSV	0.075 ppm	0.06 ppm
30 cm DSV	0.3 ppm	0.2 ppm
40 cm DSV	1.4 ppm	1.1 ppm
50x50x45 cm ³ DEV	4.0 ppm	3.1 ppm

Homogeneity - MAGNETOM Aera XQ gradients:

Homogeneity (based on highly accurate 24 plane plot)		
	Guaranteed	Typical
10 cm DSV	0.01 ppm	0.003 ppm
20 cm DSV	0.05 ppm	0.03 ppm
30 cm DSV	0.3 ppm	0.15 ppm
40 cm DSV	1.4 ppm	1.0 ppm
50x50x45 cm ³ DEV	4.0 ppm	3.0 ppm

In compliance with the German "Qualifikationsvereinbarung".

Standard deviation Vrms (Volume root-mean square) measured with highly accurate 24 plane plot method (20 points per plane).

Active shim with 3 linear and 5 non-linear channels (2nd order).

DSV = Diameter spherical volume (x, y and z direction).

DEV = Diameter elliptical volume.

Shimming	
Both: passive and active shimming. Passive shimming during installation.	
Standard active shim with 3 linear channels (1 st order). Additional 5 non linear channels (2 nd order) as option (XQ gradients).	
3D Shim	Patient-specific automated shim
	Time to shim = approx. 20 s

Shielding		
Active Shielding (AS)	5 th generation active shielding (AS) technology with counter coils	
Fringe field (axial×radial)	0.5 mT ¹	4.0 m×2.5 m
	0.1 mT	5.5 m×3.1 m
External Interference Shield (E.I.S.)	Patented shielding system integrated into the magnet	
	Continuous compensation and automatic suppression of external magnetic field interferences during measurement (caused by moving ferromagnetic objects or nearby power lines)	

Magnet cooling system	
Refill interval (typical) ¹	Not applicable
Boil-off rate (typical) ¹	0.0 l/year
Max. helium capacity	approx. 1280 liters
Minimum helium level	35%
Cryostat	Stainless steel

Hearing protection data

The A-evaluated, effective sound pressure level was measured according to NEMA MS 4-2006 (National Electrical Manufacturers Association) using the maximum gradient acoustic noise (MGAN) method.

Patient noise - MAGNETOM Aera XJ gradients:

Patient noise
Patients require hearing protection with an SNR ¹ = 18 dB or more

1. SNR = Single Number Rating

Patient noise - MAGNETOM Aera XQ gradients:

Patient noise	
Patients require hearing protection with an SNR ¹ = 22 dB or more	

1. SNR = Single Number Rating

Noise for personnel in the examination room - MAGNETOM Aera XJ gradients:

Noise for personnel in the examination room	
Noise measured	83.1 dB(A)
Hearing protection	SNR = 1 dB or more

Noise for personnel in the examination room - MAGNETOM Aera XQ gradients:

Noise for personnel in the examination room	
Noise measured	86.1 dB(A)
Hearing protection	SNR = 4 dB or more

The standard Siemens headphones are intended for communication with the patient and can be used in combination with ear plugs.

Gradient system

General features	
Actively shielded (AS) whole-body gradient coil system	
Extremely low eddy currents	
Water-cooled coil and amplifier for maximum performance	
All axes force compensated	

Gradient amplifier - MAGNETOM Aera XJ gradients:

Gradient amplifier	
- Water-cooled, highly compact, modular design - Ultra-fast solid-state technology with very low switching losses	
Max. output voltage ¹	2000 V
Max. output current ¹	625 A

Gradient amplifier - MAGNETOM Aera XQ gradients:

Gradient amplifier	
• Water-cooled, highly compact, modular design • Ultra-fast solid-state technology with very low switching losses	
Max. output voltage ¹	2250 V
Max. output current ¹	900 A

1. Values for each of the 3 gradient axes

Performance - MAGNETOM Aera XJ gradients:

Performance	
Performance for each axis	
Max. amplitude	33 mT/m
Min. rise time	264 μ s
Max. slew rate	125 T/m/s
Vector gradient performance (vector addition of all 3 gradient axes)	
Max. eff. amplitude	57 mT/m
Max. eff. slew rate	216 T/m/s
Gradient duty cycle	100%

Performance - MAGNETOM Aera XQ gradients:

Performance	
Performance for each axis	
Max. amplitude	45 mT/m
Min. rise time	225 μ s
Max. slew rate	200 T/m/s
Vector gradient performance (vector addition of all 3 gradient axes)	
Max. eff. amplitude	78 mT/m
Max. eff. slew rate	346 T/m/s
Gradient duty cycle	100%

DirectRF technology

Direct transmit technology		
Frequency stability (5 min)	$\pm 2 \times 10^{-10}$	
Frequency control	32 bits (0.015 Hz)	
Phase control	16 bits (0.006 degrees)	
Body coil	Integrated whole body no tune transmit/receive coil with 16 rungs	
	Optimized RF efficiency and signal-to-noise ratio (SNR)	
Transmitter path	Feedback loop for excellent RF stabilization	
	Transmit amplitude	16 bit control 25 ns resolution
	Gain stability (after first minute)	<0.05 dB (1 s) <0.2 dB (5 min)
Transmitter amplifier	Extremely compact, water-cooled solid state amplifier, integrated at the magnet as part of DirectRF™ technology	
	Transmit amplifier bandwidth	800 kHz
	Peak power	26.1 kW

RF receiver technology	
Number of coil elements	Up to 204
Number of independent receiver channels	24, 48, 64 ¹
Quadrature demodulation and filtering	Digital
Receiver bandwidth	500 Hz–1 MHz (for each channel)
Receiver signal resolution	32 bit
ADC sampling rate	80 MHz
Preamplifier noise figure	<0.5 dB
Dynamic range at coil connector (referred to 1 Hz resolution bandwidth)	164 dB instantaneous at receiver 169 dB with automatic gain control at local coil connector

Computer system

syngo Acquisition Workplace		
Host computer	Processor	Intel Xeon $\geq$ E5-1620 Quad-Core
	Clock rate	3.6 GHz
	Main memory (RAM)	32 GB
	1 st hard disk (system SW)	≥ 300 GB SAS
	2 nd hard disk (data base)	≥ 300 GB SAS
	3 rd hard disk (images)	≥ 300 GB SAS
	CD-R writer	Approx. 4000 images 256^2 ; DICOM Standard, ISO 9660
	DVD-R writer	Approx. 25000 images 256^2 ; DICOM Standard, ISO 9660
	Media drives	CD/DVD drive
Color LCD monitor ¹	Screen size (diagonal)	19"
	Horizontal frequency	30–100 kHz
	Vertical frequency	50–75 Hz
	Screen matrix	1280×1024 pixels
Measurement and reconstruction system (system configuration with 24 receiver channels)	Processor	Intel $\geq$ E3-1225 Quad-Core
	Clock rate	3.2 GHz, or comparable
	Main memory (RAM)	32 GB
	Hard disk for raw data	≥ 300 GB
	Hard disk for system software	≥ 100 GB
	Parallel Scan and Recon	Simultaneous scan and reconstruction of up to 12 data sets
	Reconstruction speed	12987 recons per second (256^2 FFT, full FoV) 57971 recons per second (256^2 FFT, 25% recFoV)

syngo Acquisition Workplace		
Measurement and reconstruction system (system configuration with 48 receiver channels)	Processor	Intel $\geq$ E5-2620 2x6-Core
	Clock rate	2x2.0 GHz, or comparable
	Main memory (RAM)	48 GB
	Hard disk for raw data	$\geq$ 300 GB
	Hard disk for system software	$\geq$ 100 GB
	Parallel Scan and Recon	Simultaneous scan and reconstruction of up to 12 data sets
	Reconstruction speed	22556 recons per second (256 ² FFT, full FoV) 88889 recons per second (256 ² FFT, 25% recFoV)

syngo Acquisition Workplace		
Measurement and reconstruction system (optional system configuration with 64 receiver channels)	Processor	Intel Xeon E5-2658 2x8-Core
	Clock rate	$\geq$ 2x2.1 GHz
	Main memory (RAM)	64 GB
	Hard disk for raw data	$\geq$ 400 GB
	Hard disk for system software	$\geq$ 100 GB
	Reconstruction speed	31128 recons per second (256 ² FFT, full FoV) 118519 recons per second (256 ² FFT, 25% recFoV)
	Parallel Scan and Recon	Simultaneous scan and reconstruction of up to 12 data sets
	GPGPU	1xTesla K10

syngo MR Workplace (optional)	
Color LCD monitor and host computer as for syngo Acquisition Workplace	
Euro Connector (System IEC 320)	
Load rating	100–240 Vac
	50–60 Hz
	max. 10 A

Network/data coupling

The MR system provides the level of safety according IEC 60950-1 outside patient environment. All equipment connected to the system's network/data couplings must also provide minimum level of safety according IEC 60950-1.

Cooling system

Two different customer specific cooling alternatives (Separator or Eco Chiller) are available.

Cooling system - MAGNETOM Aera XJ gradients:

Separator option (for connection to available cooling system)	Water consumption	70 l/min ¹
	Heat dissipation to water	45 kW
Eco Chiller option with automatic adaptation to the required cooling demands (e.g. different night/day mode) to decrease energy cost	GREEN Cooling Package²: • Reduces energy consumption by up to 50% • Automatic start if the surrounding temperature is 18 °C (64 °F) or less • If the temperature is less than -3 °C (27 °F) the chiller is switched off³	

Cooling system - MAGNETOM Aera XQ gradients:

Separator option (for connection to available cooling system)	Water consumption	90 l/min ¹
	Heat dissipation to water	60 kW
Eco Chiller option with automatic adaptation to the required cooling demands (e.g. different night/day mode) to decrease energy cost	GREEN Cooling Package²: • Reduces energy consumption by up to 50% • Automatic start if the surrounding temperature is 18 °C (64 °F) or less • If the temperature is less than -10 °C (14 °F) the chiller is switched off³	

1. Water temperature: 12 °C (54 °F)

2. Free Cooling Unit, optional

3. In case of clinical routine measurement conditions

Line power supply

Line power supply - MAGNETOM Aera XJ gradients:

	Values	Tolerance
Voltage	380 V, 400 V, 420 V, 440 V, 460 V, 480 V	±10%
Frequency	50/60 Hz	±1 Hz
Connection value	63 kVA for 480 V/400 V, 62 kVA for 380 V	

Line power supply - MAGNETOM Aera XQ gradients:

	Values	Tolerance
Voltage	380 V, 400 V, 420 V, 440 V, 460 V, 480 V	±10%
Frequency	50/60 Hz	±1 Hz
Connection value	88 kVA for 480 V/400 V, 90 kVA for 380 V	

For all products, line power has to be supplied via an on-site system contact or via another multipole shut-down mechanism. Room installation has to be in compliance with VDE 0100-710.

Power rating plate

Power rating plate - MAGNETOM Aera XJ gradients:

380/400/420/440/460/480	V ~
3PHASE	50/60 Hz
SHORT - TIME	kVA: 75
LONG - TIME	kVA: 20

Power rating plate - MAGNETOM Aera XQ gradients:

380/400/420/440/460/480	V ~
3PHASE	50/60 Hz
SHORT - TIME	kVA: 102
LONG - TIME	kVA: 20

Power consumption

Power consumption - MAGNETOM Aera XJ gradients:

Off ¹	6.1 kW
System ready to measure ¹	8.7 kW
Scan ¹	18.9 kW
Remarks: • All values are typical values, applicable for 400 V/50 Hz. • Consumption for optional separator pump not included.	

Power consumption - MAGNETOM Aera XQ gradients:

Off ¹	6.1 kW
System ready to measure ¹	9.3 kW
Scan ¹	22.5 kW
Remarks: • All values are typical values, applicable for 400 V/50 Hz. • Consumption for optional separator pump not included.	

1. The power consumption described herein is based on results that were achieved in a setting according to the COCIR methodology MRI - Measurement of the energy consumption (<http://www.cocir.org/site/index.php?id=46>). Since many variables impact power consumption (e.g. sequences used for scanning and sequence parameters, scan time), there can be no guarantee that each customer will achieve the same values.

Dimensions

Component	Width [cm]	Depth [cm]	Height [cm]	Weight [kg]	Heat dissipation [kW]
Examination Room					
Magnet 1.5 Tesla AS (incl. Helium)	205	137	215	3118	
Magnet in operation, incl. gradient coil, body coil, Tim table, and covers	231	405 433 ¹	219	4798	
Tim table	76	249	52–102 ² +13 mm ³		
Required min. room height clearance			240 ⁴		
Min. transport dimensions	231	155	214		
Control Room					
syngo Acquisition Workplace (table+monitor)	120	80	117 (72+45)		
Host computer	22	46	47		
syngo MR Workplace (optional) (table+monitor)	120	80	117 (72+45)		
Equipment Room					
Electronics cabinet, incl. system control, RF system, gradient power system, image processor	160	65	198 ⁵	1500	≤5 ⁶
Cooling system	65	65	189	500	

Additional system characteristics

Magnet data	
Operating field strength	1.5 Tesla
Magnet type	Superconductor
Field stability over time	<0.1 ppm/h
Weight (with cryogenics)	3118 kg
Magnet length	137 cm
System length cover to cover	145 cm
Open Bore design ¹	70 cm
Refill interval (typical) ²	Not applicable
Boil-off rate (typical) ²	0.0 l/year
Max. helium capacity	approx. 1280 liters
Minimum helium level	35%
Cryostat	Stainless steel

RF data		
Peak power of transmitter amplifier	26.1 kW	
Transmitter bandwidth	800 kHz	
Receiver bandwidth	500 Hz–1 MHz (for each channel)	
RF transmit coils		
Body coil	Max. applied RF field B1+	27 μ T
	Max. specified B1+ rms	6.2 μ T
	Distance to isocenter, where RF transmit field is reduced by	
	3 dB	0.18 m
	10 dB	0.38 m
CP Extremity	Max. applied RF field B1+	47 μ T
	Max. specified B1+ rms	17 μ T
TxRx 15 Ch Knee	Max. applied RF field B1+	47 μ T
	Max. specified B1+ rms	13 μ T
TxRx Knee 15	Max. applied RF field B1+	47 μ T
	Max. specified B1+ rms	13 μ T
TxRx CP Head	Max. applied RF field B1+	47 μ T
	Max. specified B1+ rms	12.5 μ T

Patient handling

The patient table is available in two different configurations:

Tim table		
Max. patient weight for vertical and horizontal table movement		250 kg
Max. scan range		140 cm; optional 205 cm ¹
Vertical table movement	Range	52–102 cm ² +13 mm ³
	Speed	6 cm/s
Horizontal table movement	Max. range	261 cm
	Max. speed	20 cm/s
	Horizontal accuracy for repositioning from one direction	±0.5 mm
Continuous table movement during scan capable		

Tim Dockable Table		
Max. patient weight for vertical and horizontal table movement		250 kg
Max. scan range		140 cm; optional 205 cm ¹
Vertical table movement	Range	56–106 cm ² +13 mm ³
	Speed	6 cm/s
Horizontal table movement	Max. range	261 cm
	Max. speed	20 cm/s
	Horizontal accuracy for repositioning from one direction	±0.5 mm
Continuous table movement during scan capable		

Patient comfort	
Open bore design ¹	70 cm
In-bore lighting	can be set at 6 different levels
In-bore ventilation	can be set at 6 different levels
In-bore intercom	including loudspeaker, microphone and earphones

5.4.3 syngo MR XA30A

Magnet System

Magnet parameters

Operating field strength	1.5 T
Magnet type	Superconductor
Field stability over time	< 0.1 ppm/h
Weight (with cryogenics)	3118 kg
Magnet length	1370 mm ± 3 mm
System length cover to cover	1450 mm
Open Bore design ²¹	700 mm ± 1%
Type of installation	Fixed
Decay characteristics from full field to 20 mT	approx. 20 s

Homogeneity XJ gradients(based on highly accurate 24 plane plot)

	Guaranteed	Typical
10 cm DSV	0.02 ppm	0.01 ppm
20 cm DSV	0.075 ppm	0.06 ppm
30 cm DSV	0.3 ppm	0.2 ppm
40 cm DSV	1.4 ppm	1.1 ppm
50 × 50 × 45 cm ³ DEV	4.0 ppm	3.1 ppm

In compliance with the German "Qualifikationsvereinbarung". Standard deviation Vrms (Volume root-mean square) measured with highly accurate 24 plane plot method (20 points per plane). Standard active shim with 3 linear channels (1st order). DSV = Diameter spherical volume (x, y and z direction). DEV = Diameter elliptical volume.

Homogeneity XQ gradients (based on highly accurate 24 plane plot)

	Guaranteed	Typical
10 cm DSV	0.01 ppm	0.003 ppm

Homogeneity XQ gradients (based on highly accurate 24 plane plot)

20 cm DSV	0.05 ppm	0.03 ppm
30 cm DSV	0.3 ppm	0.15 ppm
40 cm DSV	1.4 ppm	1.0 ppm
50 × 50 × 45 cm ³ DEV	4.0 ppm	3.0 ppm

In compliance with the German "Qualifikationsvereinbarung". Standard deviation Vrms (Volume root-mean square) measured with highly accurate 24 plane plot method (20 points per plane). Standard active shim with 3 linear channels (1st order) and 5 non linear channels (2nd order). DSV = Diameter spherical volume (x, y and z direction). DEV = Diameter elliptical volume.

Shimming

Both: passive and active shimming. Passive shimming during installation.

XJ and XQ:

- Standard active shim with 3 linear channels (1st order)

XQ only:

- 5 non-linear channels (2nd order)

3D Shim	Patient-specific automated shim	
	Time to shim = approx. 20 s	

Shielding

Active Shielding (AS)	5th generation active shielding (AS) technology with counter coils	
Fringe field (axial × radial)	0.5 mT ²⁾	4.0 m x 2.5 m
	0.1 mT	5.5 m x 3.1 m
External Interference Shield (E.I.S.)	Patented shielding system integrated into the magnet	
	Continuous compensation and automatic suppression of external magnetic field interferences during measurement (caused by moving ferromagnetic objects or nearby power lines)	

Magnet cooling system

Refill interval (typical) ³⁾	Not applicable
---	----------------

Magnet cooling system

Boil-off rate (typical) ⁽⁴⁾	0.0 l/year
Max. helium capacity	approx. 1280 l
Minimum helium level	450 l
Cryostat	Stainless steel

Static magnetic field

Maximum static magnetic field	1.9 T
Maximum spatial gradient of static magnetic field	11 T/m
Maximum product of the static magnetic field and spatial gradient of static magnetic field	17 T ² /m

Patient handling

Fixed patient table

Max. patient weight for vertical and horizontal table movement		250 kg (550 lbs)
Max. scan range		1400 mm optional 2050 mm ⁽⁵⁾ ±5 mm
Vertical table movement	Range	520 mm–1020 mm ⁽⁶⁾ ± 13 mm ⁽⁷⁾
Horizontal table movement	Max. range	2619 mm ±10 mm
	Max. speed	200 mm/s ±2 mm/s
	Horizontal accuracy for repositioning from one direction	± 0.5 mm

Continuous table movement during scan capable

Dockable patient table

Max. patient weight for vertical and horizontal table movement		250 kg (550 lbs)
Max. scan range		1400 mm optional 2050 mm ⁸⁾ ±5 mm
Vertical table movement	Range	560 mm–1060 mm ⁹⁾ ± 13 mm ¹⁰⁾
Horizontal table movement	Max. range	2619 mm ±10 mm
	Max. speed	200 mm/s ±2 mm/s
	Horizontal accuracy for repositioning from one direction	± 0.5 mm
Continuous table movement during scan capable		

Hearing protection data

The A-evaluated, effective sound pressure level was measured according to NEMA MS 4-2010 (National Electrical Manufacturers Association) using the maximum gradient acoustic noise (MGAN) method.

Patient noise

XJ gradients	Patients require hearing protection with an SNR = 18 dB or more ¹¹⁾
XQ gradients	Patients require hearing protection with an SNR = 22 dB or more

Noise for personnel in the examination room

XJ gradients	Noise measured	83.1 dB(A)
	Hearing protection	SNR = 1 dB or more
XQ gradients	Noise measured	86.1 dB(A)
	Hearing protection	SNR = 4 dB or more

Gradient system

- Actively shielded (AS) whole-body gradient coil system
- Extremely low eddy currents
- Water-cooled coil and amplifier for maximum performance
- All axes force compensated

XJ gradients

Gradient amplifier¹²⁾

- Water-cooled, highly compact, modular design
- Ultra-fast solid-state technology with very low switching losses

Max. output voltage	2000 V
Max. output current	625 A
Max. power	1.250 MW

XQ gradients

Gradient amplifier¹³⁾

- Water-cooled, highly compact, modular design
- Ultra-fast solid-state technology with very low switching losses

Max. output voltage	2250 V
Max. output current	900 A
Max. power	2.025 MW

XJ gradients

Gradient performance for each axis

Max. amplitude	33 mT/m
Min. rise time	264 µs
Max. slew rate	125 T/m/s

Vector gradient performance (vector addition of all 3 gradient axes)

Max. eff. amplitude	57 mT/m
Max. eff. slew rate	216 T/m/s
Gradient duty cycle	100%

XQ gradients

Gradient performance for each axis

Max. amplitude	45 mT/m
Min. rise time	225 μ s
Max. slew rate	200 T/m/s

Vector gradient performance (vector addition of all 3 gradient axes)

Max. eff. amplitude	78 mT/m
Max. eff. slew rate	346 T/m/s
Gradient duty cycle	100%

RF technology

Transmit technology

Frequency stability (5 min)	$\pm 2 \times 10^{-10}$	
Frequency control	32 bits (0.03 Hz)	
Phase control	16 bits (0.006°)	
Body coil	• Integrated whole body no tune transmit/receive coil with 16 rungs • Optimized RF efficiency and signal-to-noise ratio (SNR) • Real-time Feedback loop for unmatched RF stability	
Transmitter path	Feedback loop for excellent RF stabilization	
	Transmit amplitude	16 bit control 25 ns resolution
	Gain stability (after first minute)	< 0.05 dB (1 s) < 0.2 dB (5 min)
Transmit amplifier	Extremely compact, water-cooled solid state amplifier, fully integrated at the magnet	
	Transmit amplifier bandwidth	500 kHz
Peak power of transmitter amplifier	26.1 kW	

RF receiver technology

Maximum number of channels ¹⁴⁾	204
Number of independent receiver channels that can be used simultaneously in one single scan in one FOV, each generating an independent partial image	24 48 64 ¹⁵⁾
Quadrature demodulation and filtering	Digital
Receiver bandwidth	500 Hz–1 MHz (for each channel)
Receiver signal resolution	32 bit
ADC sampling rate	80 MHz
Preamplifier noise figure	< 0.5 dB
Dynamic range at coil connector (referred to 1 Hz resolution bandwidth)	164 dB instantaneous at receiver 169 dB with automatic gain control at local coil connector

Computer system

The specifications contained in this section apply to new system deliveries with XA30A and depending on appropriate hardware preconditions also to systems which are upgraded to XA30A.

Host computer

Processor	Intel Xeon ≥ W-2133 (6 Core)
Clock rate	3.6 GHz, or comparable
Main memory (RAM)	64 GB
Hard disk (DICOM Standard, ISO 9660)	SSD: 480 GB Approx. images: 256 ²²⁰⁾
CD/DVD drive	Not built in, but optionally connectable by USB
Media drive	SDHC card reader

Advanced Host computer

Processor	1×Intel Xeon ≥ W-2145 (8 Core)
Clock rate	3.7 GHz, or comparable
Main memory (RAM)	96 GB
Hard disk (DICOM Standard, ISO 9660)	SSD: 480 GB Approx. images: 256 ²²¹⁾
CD/DVD drive	Not built in, but optionally connectable by USB
Media drive	SDHC card reader

Color LCD monitor²²⁾

High resolution widescreen monitor:²³⁾

- Horizontally tiltable, forward and backward
- Automatic backlight control for long-term brightness stability

Screen size (diagonal)	24"
Horizontal scanning frequency	31 kHz–76 kHz
Vertical scanning frequency	59 Hz–61 Hz
Screen matrix	1920 pixels × 1200 pixels

syngo MR Workplace

Processor	Intel Pentium J5005 (4 Core)
Clock rate	1.5 GHz
Main memory (RAM)	8 GB
Hard disk	SSD: 128 GB
CD/DVD drive	Not built in, but optionally connectable by USB
Audio	On board
Ethernet	On board, 1 Gbit

Measurement and reconstruction system for Tim [204 × 24]

Processor	1 × Intel Xeon E3-1275v5 (4 Core)
Clock rate	3.6 GHz, or comparable
Main memory (RAM)	48 GB
SSD for raw data	≥ 480 GB
SSD for system software	≥ 240 GB
Reconstruction speed	16 949 recons per second (256 ² FFT, full FoV) 78 431 recons per second (256 ² FFT, 25% recFoV)
Parallel Scan and Recon	Simultaneous scan and reconstruction of up to 12 data sets

Advanced measurement and reconstruction system for Tim [204 × 24]

Processor	2 × Intel Xeon E5-2620v4 (8 Core)
Clock rate	2 × 2.1 GHz, or comparable
Main memory (RAM)	96 GB
SSD for raw data	≥ 480 GB
SSD for system software	≥ 240 GB
Reconstruction speed	40 404 recons per second (256 ² FFT, full FoV) 149 532 recons per second (256 ² FFT, 25% recFoV)
Parallel Scan and Recon	Simultaneous scan and reconstruction of up to 12 data sets
GPGPU	1 × Nvidia Quadro P2000

Measurement and reconstruction system for Tim [204 × 48]

Processor	2 × Intel Xeon E5-2609v4 (8 Core)
Clock rate	2 × 1.7 GHz, or comparable
Main memory (RAM)	48 GB
SSD for raw data	≥ 480 GB
SSD for system software	≥ 240 GB
Reconstruction speed	32 854 recons per second (256 ² FFT, full FoV) 122 137 recons per second (256 ² FFT, 25% recFoV)
Parallel Scan and Recon	Simultaneous scan and reconstruction of up to 12 data sets

Advanced measurement and reconstruction system for Tim [204 × 48]

Processor	2 × Intel Xeon E5-2620v4 (8 Core)
Clock rate	2 × 2.1 GHz, or comparable
Main memory (RAM)	96 GB
SSD for raw data	≥ 480 GB
SSD for system software	≥ 240 GB
Reconstruction speed	40 404 recons per second (256 ² FFT, full FoV) 149 532 recons per second (256 ² FFT, 25% recFoV)
Parallel Scan and Recon	Simultaneous scan and reconstruction of up to 12 data sets
GPGPU	1 × Nvidia Quadro P2000

Measurement and reconstruction system for Tim [204 × 64]

Processor	2 × Intel Xeon E5-2620v4 (8 Core)
Clock rate	2 × 2.1 GHz, or comparable
Main memory (RAM)	64 GB
SSD for raw data	≥ 480 GB
SSD for system software	≥ 240 GB
Reconstruction speed	40 404 recons per second (256 ² FFT, full FoV) 149 532 recons per second (256 ² FFT, 25% recFoV)
Parallel Scan and Recon	Simultaneous scan and reconstruction of up to 12 data sets
GPGPU	1 × Nvidia Quadro P4000

Advanced measurement and reconstruction system for Tim [204 × 64]

Processor	2 × Intel Xeon E5-2640v4 (10 Core)
Clock rate	2 × 2.4 GHz, or comparable
Main memory (RAM)	128 GB
SSD for raw data	≥ 2 × ≥ 480 GB
SSD for system software	≥ 240 GB
Reconstruction speed	52 493 recons per second (256 ² FFT, full FoV) 224 719 recons per second (256 ² FFT, 25% recFoV)
Parallel Scan and Recon	Simultaneous scan and reconstruction of up to 12 data sets
GPGPU	2 × Nvidia Quadro P4000

Euro Connector (System IEC 320)

Load rating	100 Vac–240 Vac
	50 Hz–60 Hz
	max. 10 A

Cooling system

Two different customer specific cooling alternatives (Separator or Eco Chiller) are available.

Separator option (for connection to available cooling system)	Water consumption	XJ gradients: 80 l/min ± 10 l/min ²⁹⁾ XQ gradients: 100 l/min ± 10 l/min ²⁹⁾
	Heat dissipation to water	XJ gradients: 45 kW XQ gradients: 60 kW
Eco Chiller option with automatic adaptation to the required cooling demands (e.g. different night/day mode) to decrease energy cost ²⁹⁾	GREEN Cooling Package: ³⁰⁾ - Automatic start if the surrounding temperature is 18°C (64°F) or less - If the temperature is less than -10°C (14°F) the chiller is switched off ³¹⁾	

Line power supply

XJ gradients

Line power supply

	Values	Tolerance
Voltage	380 V, 400 V, 420 V, 440 V, 460 V, 480 V	± 10%
Frequency	50 Hz/60 Hz	± 1 Hz
Connection value	63 kVA	

XQ gradients

Line power supply

	Values	Tolerance
Voltage	380 V, 400 V, 420 V, 440 V, 460 V, 480 V	± 10%
Frequency	50 Hz/60 Hz	± 1 Hz
Connection value	88 kVA	

For all products, line power has to be supplied via an on-site system contact or via another multipole shut-down mechanism. Room installation has to be in compliance with VDE 0100-710.

In all countries, compliance with local and national legal regulations is required. However, we strongly recommend compliance with the regulations described herein – to the extent permitted by relevant local and national laws – in order to ensure the safety of operating personnel, patients, and third parties.

Power rating plates

XJ gradients

380/400/420/440/460/480 V ~
3PHASE 50/60 Hz
SHORT – TIME kVA: 75
LONG – TIME kVA: 20

XQ gradients

380/400/420/440/460/480 V ~
3PHASE 50/60 Hz
SHORT – TIME kVA: 102
LONG – TIME kVA: 20

Power consumption

XJ gradients

Power consumption³²⁾

System off	6.1 kW
System ready to measure	8.7 kW
Scan	18.9 kW

XQ gradients

Power consumption³³⁾

System off	6.1 kW
System ready to measure	9.3 kW
Scan	22.5 kW

Dimensions

Examination Room

Component	Width	Depth	Height	Weight
Magnet 3 T AS (incl. Helium)	2.05 m	1.37 m	2.15 m	3118 kg
Magnet in operation, incl. gradient coil, body coil, fixed patient table (Tim table), and covers	2.31 m	4.05 m 4.33 m ³⁴⁾	2.19 m	4798 kg
Fixed patient table (Tim table)	0.76 m	2.49 m	0.52 m–1.02 m; ³⁵⁾ + 0.013 m ³⁶⁾	

Examination Room

Component	Width	Depth	Height	Weight
Required min. room height clearance			2.4 m ³⁷⁾	
Min. transport dimensions	2.31 m	1.55 m	2.14 m	

Control Room

Component	Width	Depth	Height	Weight
syngo Acquisition Workplace (table+monitor)	1.2 m	0.8 m	1.17 m (0.72 m + 0.45 m)	
Host computer	0.22 m	0.46 m	0.47 m	
syngo MR Workplace (optional) (table+monitor)	1.2 m	0.8 m	1.17 m (0.72 m + 0.45 m)	

Equipment Room

Component	Width	Depth	Height	Weight
Electronics cabinet, incl. system control, RF system, gradient power system, image processor	1.6 m	0.65 m	1.98 m ³⁸⁾	1500 kg
Cooling system	0.65 m	0.65 m	1.89 m	500 kg
Heat dissipation	≤ 5 kW ³⁹⁾			

Additional system characteristics

Transmit technology

Frequency stability (5 min)	$\pm 2 \times 10^{-10}$	
Frequency control	32 bits (0.03 Hz)	
Phase control	16 bits (0.006°)	
Body coil	• Integrated whole body no tune transmit/receive coil with 16 rungs • Optimized RF efficiency and signal-to-noise ratio (SNR) • Real-time Feedback loop for unmatched RF stability	
Transmitter path	Feedback loop for excellent RF stabilization	
	Transmit amplitude	16 bit control 25 ns resolution
	Gain stability (after first minute)	$< 0.05 \text{ dB (1 s)}$ $< 0.2 \text{ dB (5 min)}$
Transmit amplifier	Extremely compact, water-cooled solid state amplifier, fully integrated at the magnet	
	Transmit amplifier bandwidth	500 kHz

Transmit technology

Peak power of transmitter amplifier	26.1 kW
-------------------------------------	---------

RF receiver technology

Maximum number of channels ⁴⁾	204
Number of independent receiver channels that can be used simultaneously in one single scan in one FOV, each generating an independent partial image	24 48 64 ⁵⁾
Quadrature demodulation and filtering	Digital
Receiver bandwidth	500 Hz–1 MHz (for each channel)
Receiver signal resolution	32 bit
ADC sampling rate	80 MHz
Preamplifier noise figure	< 0.5 dB
Dynamic range at coil connector (referred to 1 Hz resolution bandwidth)	164 dB instantaneous at receiver 169 dB with automatic gain control at local coil connector

RF transmit coils

Body coil	Max. applied RF field B1+	27 μ T
	Max. specified B1+ rms	6.2 μ T
	Distance to isocenter, where RF transmit field is reduced by	
	3 dB	0.18 m
	10 dB	0.38 m
TxRx CP Extremity	Max. applied RF field B1+	47 μ T
	Max. specified B1+ rms	17 μ T
TxRx 15 Ch Knee	Max. applied RF field B1+	47 μ T
	Max. specified B1+ rms	13 μ T

RF transmit coils

TxRx Knee 15	Max. applied RF field B1+	47 μ T
	Max. specified B1+ rms	12 μ T
TxRx Knee 18	Max. applied RF field B1+	47 μ T
	Max. specified B1+ rms	7.6 μ T
TxRx CP Head	Max. applied RF field B1+	47 μ T
	Max. specified B1+ rms	10.6 μ T

Patient handling

Fixed patient table

Max. patient weight for vertical and horizontal table movement		250 kg (550 lbs)
Max. scan range		1400 mm optional 2050 mm ⁶⁾ ± 5 mm
Vertical table movement	Range	520 mm–1020 mm ⁷⁾ ± 13 mm ⁸⁾
Horizontal table movement	Max. range	2619 mm ± 10 mm
	Max. speed	200 mm/s ± 2 mm/s
	Horizontal accuracy for repositioning from one direction	± 0.5 mm

Continuous table movement during scan capable

Dockable patient table

Max. patient weight for vertical and horizontal table movement		250 kg (550 lbs)
Max. scan range		1400 mm optional 2050 mm ⁹⁾ ± 5 mm

Vertical table movement	Range	560 mm–1060 mm ¹⁰⁾ ± 13 mm ¹¹⁾
Horizontal table movement	Max. range	2619 mm ±10 mm
	Max. speed	200 mm/s ±2 mm/s
	Horizontal accuracy for repositioning from one direction	± 0.5 mm

Continuous table movement during scan capable

Patient comfort

Open bore design ¹²⁾	700 mm ± 1%
In-bore lighting	can be set at 6 different levels
In-bore ventilation	can be set at 6 different levels
In-bore intercom	including loudspeaker, microphone and earphones

5.5 Capabilities of the Medical device and ways of integration with other MDs

The following medical devices may be integrated with MAGNETOM Aera:

Medical Device	Legal Manufacturer	Way of integration
Patient table	TRUMPF Medizin Systeme GmbH + Co.KG Carl-Zeiss-Str. 7-9 07318 Saalfeld, Germany www.trumpf-med.com	Component
Respiration cushion, respiratory belt	Pi-Products GmbH Heinrich-Hertz-Straße 8a 92224 Amberg, Germany www.pi-products.de	Accessories
Application cushion for PERU/PMU	Polyform GmbH & Co. KG Kunststofftechnik Braasstraße 15 DE-31737 Rinteln, Germany www.polyform.de	Accessories


Medical Device	Legal Manufacturer	Way of integration
Shoulder Small/ Large 16, Shoulder Shape 16, Hand/Wrist 16, Foot/Ankle 16 Special Purpose 4 Head/Neck 16 Spine 24, Body 6 TxRx Knee 18, Flex Coil Interface, Tim Coil Interface	Siemens Shenzhen Magnetic Resonance Ltd. (SSMR) Siemens MRI Center Gaoxin C. Ave., 2nd Hi-Tech Industrial Park Shenzhen 518057, P.R. China www.siemens.com/healthineers Authorized Representative (MDD 93/42/EEC) Siemens Healthcare GmbH Henkestr. 127 91052 Erlangen, Germany	Coil
Prostate 2, 2-/4-/8-Channel Sentinelle Breast Coil, 2-/10-/16-Channel Sentinelle Breast Coil	Invivo Corporation 3545 SW 47th Ave Gainesville, Florida 32608, USA Authorized Representative (MDD93/42/EEC) Philips Medical Systems Nederland B.V. Veenpluis 6 5684 PC, The Netherlands	Coil
TxRx 15 Ch Knee/ TxRx CP Extremity TxRx Knee 15 TxRx CP Head Pediatric 16 TxRx Knee 15 Flare Contour 24 coil	Quality Electrodynamics, LLC 6655 Beta Drive Suite 100, Mayfield Village, OH 44143, USA https://qedinnovations.com Authorized Representative (MDD 93/42/EEC) EMERGO EUROPE Prinsessegracht 20 2514 AP The Hague, Netherlands	Coil
Endorectal	Bayer Medical Care, Inc. 1 Bayer Drive Indianola, PA 15051, USA Authorized Representative (MDD 93/42/EEC) Bayer Medical Care BV Horsterweg 24 6199 AC Maastricht Airport, Netherlands http://pharma.bayer.com	Coil

Medical Device	Legal Manufacturer	Way of integration
4 Ch BI Breast Coil Breast BI 7	NORAS MRI products GmbH Leibnizstrasse 4 D-97204 Höchberg, Germany www.noras.de	Coil
Positioning aids/ Table pads	Polyform GmbH & Co. KG Kunststofftechnik Braasstraße 15 DE-31737 Rinteln, Germany www.polyform.de	Accessories
Sandbags, RF homogeneity cushion, vacuum cushion, belts, positioning aids for Breast 18	Pi-Products GmbH Heinrich-Hertz Straße 8a 92224 Amberg, Germany www.pi-products.de	Accessories

5.6 Labeling information and requirements


The following labels can be found at the electronics cabinet:

Product identification label (MDD)

MDD = Medical Devices Directive

Product identification label (MDR)

MDR = Medical Device Regulation

(1) Configurable device


A configurable device is a device that consists of several components which can be assembled by the manufacturer in multiple configurations. Those individual components may be devices in themselves.

(2) Identification as medical device

(3) Serial number of the system

(4) Unique device identification

Combination label (MDR)

Labels according to EN 50419:2006

Within the EU, products identified with this symbol are subject to guideline 2002/96/EC for old electrical or electronics system and 2012/19/EU respectively. Please contact Siemens Service in case of questions about returning and disposing the MR system and/or its components and accessories

Approval identification for Canada/USA

CSA = Canadian Standards Association

SGS North America

Revision label

Revision:														Model - No.:																		
X	X	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65
66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98

The current revision is marked with an X in the respective field. The example above shows revision level 01 of the respective component.


The following labels can be found on the magnet cover:

Product identification label (MDD)

MDD = Medical Devices Directive

Product identification label (MDR)

MDR = Medical Device Regulation

(1) Configurable device


A configurable device is a device that consists of several components which can be assembled by the manufacturer in multiple configurations. Those individual components may be devices in themselves.

(2) Identification as medical device

(3) Serial number of the system

(4) Unique device identification

Combination Label (MDR)

Approval identification for Canada/USA

CSA = Canadian Standards Association

SGS North America

The following labels can be found on the magnet bore (front)

Risk of injury

Observe operator manual


For Laser light localizer:


Safety certificate label (International)



Safety certificate label (USA only)




The following labels can be found on the coils:


Coil Name	Labels
Head/Neck coils	Please note: Siemens AG until mid of July 2016 (also for refurbished systems) Siemens Healthcare GmbH starting mid of July 2016 Head/Neck 16 MR Coil   Head/Neck 20 MR Coil  
Spine coils	Spine 24 MR Coil 1.5T   Spine 32 MR Coil 1.5T 

	
Tx/Rx CP Head Coil for vascular examinations	 
Body coils	Body 6 MR Coil 1.5T   Body 6 long MR Coil 1.5T 

	Body 18 MR Coil 1.5T   Body 18 long MR Coil 1.5T  
Body Coils for expert examinations	Body 30 MR Coil 1.5T 

	 ES: 01 02 03 04 05 06 07 08 09 10 11 12 13 14 15 16 17 18 19 SIEMENS (TP) Model No. 10835444 (S) Serial No. 1001 IVK Siemens AG, Wittenbergplatz 2, DE - 80333 München Made in Germany 5004313
	Contour 24 coil  Quality Electrodynamics, LLC 6055 Beta Drive, Suite 100, Mayfield Village, OH 44130 USA Distributed by Siemens Healthcare GmbH Contour 24 (Receive Only) for Siemens 3.0T, 63.6 MHz Systems Siemens (MFP) (TP) 11134546 QED (MFP) 07000283 QED (S) 140104 ES: 01 02 03 04 05 06 07 08 09 10 11 12 13 14 15 16 17 18 19 CE 0086 MCSG GmbH Schiffgraben 41 30175 Hannover Germany IVK GTIN CTN PROD DATE SERIAL MODEL ORIGIN Conforms to US Std 40843.1, ANSI/AAMI Std ES60601-1 IEC Std 60601-1, 60601-2-33 Certificate: CDR, Std C2.2, No. 602.1.8.60601.1 Label No. 6000052
Peripheral Angio 36 Coil	 SIEMENS PA 36 MR Coil 1.5T (R) 8405080811495 GUN (TH) 1049542 MODEL (AZ) 276 ORIGIN DE (S) 1432 SERIAL IVK Siemens Healthcare GmbH, Munchen 127 80333 München, Germany Made in Germany  SIEMENS (TP) Model No. 1049542 (S) Serial No. 1432 IVK Siemens AG, Wittenbergplatz 2, DE - 80333 München Made in Germany ES: 01 02 03 04 05 06 07 08 09 10 11 12 13 14 15 16 17 18 19
Tx/Rx Knee 15 coils for joint examinations	Tx/Rx Knee 15 Flare Coil  Quality Electrodynamics, LLC 6055 Beta Drive, Suite 100 Mayfield Village, OH 44130 USA Distributed by Siemens Healthcare GmbH Tx/Rx Knee 15 Flare MR Coil 1.5T Siemens (MFP) (TP) 10675089 QED (MFP) 07000147 QED (S) 140104 ES: 01 02 03 04 05 06 07 08 09 10 11 12 13 14 15 16 17 18 19 CE 0086 MCSG GmbH Schiffgraben 41 30175 Hannover Germany IVK GTIN CTN PROD DATE SERIAL MODEL ORIGIN Conforms to US Std 40843.1, ANSI/AAMI Std ES60601-1 IEC Std 60601-1, 60601-2-33 Certificate: CDR, Std C2.2, No. 602.1.8.60601.1 Label Part No. 6000740

	Tx/Rx Knee 15 Coil
Tx/Rx Knee 18 Coil	TX/RX Knee 18
Tx/Rx CP Extremity Coil	SIEMENS CP Extremity MR Coil SIEMENS
Foot/Ankle Coil	000 1.5 T Foot Ankle Coil

Hand/Wrist Coil	 
Shoulder Coils	Shoulder 16 Large   Shoulder 16 Small  
Shoulder Shape 16 Coil	
Breast 18 Coil for expert examinations	

Breast Biopsy Coils 4-channel BI Breast Coil	4 Ch BI Breast MR Coil KR 1.5T Modular equipment for MRI imaging Manufacturer: Noras MRI products GmbH Leibnizstr. 4 97204 Höchberg Germany Distributor: Siemens Healthcare GmbH Model No.: 110121 Rev. 07 CE 0123 4 Ch BI Breast MR Coil KR 1.5T / 63,6 MHz Manufacturer: Noras MRI products GmbH Leibnizstr. 4 97204 Höchberg Germany Distributor: Siemens AG Wilhelmsbader Platz 2 80333 München Germany Part No.: 10195524 Serial No.: 1001 Plate No.: 05 Quantity: 05 CE 0123 Insertion Plate Breast Coil 1.5T Modular equipment for MRI imaging Manufacturer: Noras MRI products GmbH Leibnizstr. 4 97204 Höchberg Germany Distributor: Siemens Healthcare GmbH Model No.: 110121 Rev. 06 CE 0123 Insertion Plate Breast Coil 1.5T / 63,6 MHz Manufacturer: Noras MRI products GmbH Leibnizstr. 4 97204 Höchberg Germany Distributor: Siemens AG Wilhelmsbader Platz 2 80333 München Germany CE 0123 Patient Pad Breast Coil 1.5T Modular equipment for MRI imaging Manufacturer: Noras MRI products GmbH Leibnizstr. 4 97204 Höchberg Germany Distributor: Siemens Healthcare GmbH Model No.: 110121 Rev. 07 CE 0123 Patient Pad Breast Coil 1.5T / 63,6 MHz Manufacturer: Noras MRI products GmbH Leibnizstr. 4 97204 Höchberg Germany Distributor: Siemens AG Wilhelmsbader Platz 2 80333 München Germany CE 0123
Modular Breast Coils for visualization and biopsy 2-,4-, 8 ch Sentinelle Coil 2-,10-, 16 ch Sentinelle Coil	Sentinelle by InVivo 2/4CH BREAST SEN MR COIL KIT 1.5T REF (1P) 10 60 66 02 Distributed by: Siemens Healthcare GmbH Henkestr. 127 91052 Erlangen Germany Rev (2P) 03

Sentinelle
by Invivo

CE 0123

2/4CH BREAST SEN MR COIL KIT 1.5T

Distributed by:
Siemens AG
Wittelsbacherplatz 2
DE-80333 Muenchen
Germany

REF (1P) 10 60 66 02

SN (S) XXXXXX

Rev (2P) 02

GBR Authorized Representative Ltd
The White House
2 Watlington Road
Goslington, Northampton NN4 6JF
United Kingdom
+44 (0)1933 555 721
www.invivo-gbr.com

Sentinelle
by Invivo

CE 0123

8CH BREAST SEN MR COIL KIT 1.5T

Distributed by:
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany

REF (1P) 10 60 66 01

SN (S) serial

Rev (2P) 06

Invivo Corporation
3045 304 Avenue
Beverly Hills, CA 90210-1699
U.S.A. +1 (310) 355-8810

Sentinelle
by Invivo

CE 0123

8CH BREAST SEN MR COIL KIT 1.5T

Distributed by:
Siemens AG
Wittelsbacherplatz 2
DE-80333 Muenchen
Germany

REF (1P) 10 60 66 01

SN (S) XXXXXX

Rev (2P) 06

GBR Authorized Representative Ltd
The White House
2 Watlington Road
Goslington, Northampton NN4 6JF
United Kingdom
+44 (0)1933 555 721
www.invivo-gbr.com

Sentinelle
by Invivo

CE 0123

16 CH BREAST SEN MR COIL KIT 1.5T

Distributed by:
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany

REF (1P) 10 60 66 09

SN (S) serial

Rev (2P) 03

GBR Authorized Representative Ltd
The White House
2 Watlington Road
Goslington, Northampton NN4 6JF
United Kingdom
+44 (0)1933 555 721
www.invivo-gbr.com

Sentinelle
by Invivo

CE 0123

16 CH BREAST SEN MR COIL KIT 1.5T


Distributed by:
Siemens AG
Wittelsbacherplatz 2
DE-80333 Muenchen
Germany


REF (1P) 10 60 66 09

SN (S) serial


Rev (2P) 03

GBR Authorized Representative Ltd
The White House
2 Watlington Road
Goslington, Northampton NN4 6JF
United Kingdom
+44 (0)1933 555 721
www.invivo-gbr.com

Breast Coil for visualization and biopsy Breast BI7 MR Coil	 Breast BI 7 MR Coil 1.5T GTIN Origin Rev Serial Model REF 118879 Rev. 01 SN XXXX 2019-07-08 CE 0123
Prostate 2	 Sentinelle by InVivo PROSTATE 2 SEN MR COIL KIT 1.5T REF (1P) 10 60 66 24 SN (S) serial Rev (2P) 002 CE 0123 Distributed by: Siemens Healthcare GmbH Henkelstr. 127 91052 Erlangen Germany CEM Authorized Representative Ltd The White House 2 Meadway Cuddingbury, OLT 396 United Kingdom +44 (0)1775 539 721 www.mfa-gbm.com
	 Sentinelle by InVivo PROSTATE 2 SEN MR COIL KIT 1.5T REF (1P) 10 60 66 24 SN (S) XXXXXX Rev (2P) 002 CE 0123 Distributed by: Siemens AG Waldschloßplatz 2 DE-80333 München Germany CEM Authorized Representative Ltd The White House 2 Meadway Cuddingbury, OLT 396 United Kingdom +44 (0)1775 539 721 www.mfa-gbm.com For Rx Use Only
Special Purpose Coils	 4-Channel Special - Purpose Coil (01) 0405000014808 GTIN (140) 10499077 MODEL (422) 156 ORIGIN CN (21) 1025 SERIAL IVK Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center, Gaoxin C. Ave., 2nd, 10 - Tech Industrial Park, 518057 Shenzhen, China Made in China
	 SIEMENS (1P) MODEL No. 10499077 (2S) SERIALLOT No. 1017 MADE IN CHINA IVK PS 02 03 04 05 06 07 08 09 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24 25 26 27 28 MR2 AAN330071 CE 0123
Flex Large 4 Coil	 SIEMENS Flex Large 4 MR coil 1.5T (01) 0405000014784 GTIN (140) 10105532 MODEL (422) 276 ORIGIN DE (21) 30036 SERIAL IVK Siemens Healthcare GmbH, Henkelstr. 127 91052 Erlangen, Germany Made in Germany

	
Flex Small 4 Coil	
	
UltraFlex Large 18 Coil	
UltraFlex Small 18 Coil	
Small Loop Coil	
	

Medium Loop Coil	 
Large Loop Coil	 
Prostate Endo Coil	
Colorectal Endo Coil	
Cervix Endo Coil	

	 medrad® eCoil MR Endorectal Coil Antenne IRM Endorectale MEDRAD eCoil Disposable Jetable pour eCoil de MEDRAD Siayer Medical Care Inc. 1 Siayer Drive Iselin, NJ 07034-0788 USA Siayer Medical Care B.V. Rustenburg 24 6750 AC, Mook-Weert, Arnhem The Netherlands 日本メディカル株式会社 〒330-0001 大塚市北沢 4-9 日本 REF BCR15 PN 40774403 LOT 10004508 EXP 06/08 CE 0086						
Pediatric 16 Coil	 Quality Electrodynamics, LLC 6655 Beta Drive, Suite 100 Mayfield Village, OH 44143 USA Distributed by Siemens Healthcare GmbH Pediatric 16, A 1.5T T1M Coil (Receiver Only) <table border="1"> <tr> <td>Siemens REF (1P) 10835522</td> <td>03 04 05</td> </tr> <tr> <td>QED REF Q7000061</td> <td>06 07 08</td> </tr> <tr> <td>SN (S) NNNN</td> <td>09 10 11</td> </tr> </table> Freq MHz: 63.6 CE 0086 MESS GmbH Schiffgraben 41 30175 Hannover Germany IVK Rx ONLY (01)0081484920416 GTIN (11)YYMMDD PROD DATE (21)NNNN SERIAL (240)10835522 MODEL (422)840 ORIGIN Label Part No. 6002578	Siemens REF (1P) 10835522	03 04 05	QED REF Q7000061	06 07 08	SN (S) NNNN	09 10 11
Siemens REF (1P) 10835522	03 04 05						
QED REF Q7000061	06 07 08						
SN (S) NNNN	09 10 11						
	 Quality Electrodynamics, LLC 700 Beta Drive, Suite 100 Mayfield Village, OH 44143 USA Distributed by SIEMENS Pediatric 16, A 1.5T T1M Coil (Rx Only) <table border="1"> <tr> <td>Siemens REF (1P) 10835522</td> <td>03 04 05</td> </tr> <tr> <td>QED REF Q7000061</td> <td>06 07 08</td> </tr> <tr> <td>SN (S) 1001</td> <td>09 10 11</td> </tr> </table> Freq MHz: 63.6 CE 0086 MESS GmbH Schiffgraben 41 30175 Hannover Germany IVK Rx ONLY (01)0081484920416 GTIN (11)YYMMDD PROD DATE (21)NNNN SERIAL (240)10835522 MODEL (422)840 ORIGIN Label Part No. 6002578	Siemens REF (1P) 10835522	03 04 05	QED REF Q7000061	06 07 08	SN (S) 1001	09 10 11
Siemens REF (1P) 10835522	03 04 05						
QED REF Q7000061	06 07 08						
SN (S) 1001	09 10 11						

* Siemens AG (as the Manufacturer) may be indicated on the original labeling of the Manufacturer for the systems which were produced before middle of July 2016 (also for refurbished systems);

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~3/PE 380/400/420/440/460/480 В

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GTIN	Global number GTIN
Model	Model
	Serial number
	Manufacturer
	Catalog number
	Not to throw out. To give to utilization point.
	Conformity of EC
	EC Representative
	Caution! Refer to instruction for use
	Refer to instruction for use
	Class II
	Applied part type BF

6. Safety Requirements (Conditions and rules of operation)

6.1 Terms and Definitions

Term used	Explanation
User/ Operator/ Operating personnel	Person who operates the system or software, takes care of the patient or reads images Typically physicians, trained radiological technicians, or technologists
System owner	Person who is responsible for the MR environment. This includes legal requirements, emergency plans, employee information and qualifications, as well as maintenance/repair.
MR worker	Person who works within the controlled access area or MR environment User/Operator as well as further personnel (for example, cleaning staff, facility manager, service personnel)
Siemens Service/service personnel	Group of specially trained persons who are authorized by Siemens to perform certain maintenance activities References to "Siemens Service" include service personnel authorized by Siemens.

6.2 Preamble

This section contains the most important safety aspects for which the owner of the MR system, will be responsible. These include legal requirements, emergency plans, employee information and qualifications, as well as requirements that must be met in the examination room.

6.3 Common causes for accidents

Of particular importance is the obligation to inform employees and contractors. Operating personnel as well as personnel who are not regulars in the examination room (e.g. cleaning personnel, rescue personnel) must be informed about the special conditions in the examination room. Magnetizable devices (e.g. floor polishers, vacuum cleaners, wheelchairs, metal gurneys) must not be used in the examination room.

Please note all safety instructions applicable to all users which are described in the operator manual of the MR system.

6.4 Responsibility

Siemens accepts no responsibility for the safety, reliability, and performance of the MR system, if the MR system is not used in accordance with the instructions for use (Operator Manual, System Owner Manual). Siemens is also not responsible for any direct or indirect damages caused by incorrect operation.

This includes, but is not limited to, accidents with ferromagnetic objects. This applies even if the consequences only become obvious at a later point in time.

6.5 Legal regulations

6.5.1 Country-specific regulations

Local and national legal regulations must be observed. It is the operator's responsibility to follow local statutory requirements regarding access to the controlled access area.

The local regulations also define the acceptable exposure limits regarding noise as well as magnetic stray fields for users and patients.

6.5.2 Pressure Equipment Directive

The super-conductive magnet is classified as pressure equipment. National guidelines for starting up and operating pressure equipment must be observed.

In Europe, the Pressure Equipment Directive (97/23/EG) regulates the sale of pressure devices.

In Germany both the Pressure Equipment Directive as well as the relevant Occupational Safety Regulations (BetrSichV) for system start-up and operation apply.

6.5.3 Electromagnetic fields

The 0.5 mT line in the examination room defines the controlled access area of the static magnetic field. For controlling access to this area adequate rules must be established. The potential risks from the attraction of magnetizable objects or from torque on such materials must also be considered. The possibility must be addressed that persons with medical implants, such as pacemakers, could enter the area without realizing that the medical implants may be adversely affected by the field.

Outside the controlled access area, electromagnetic interferences meet the requirements according to IEC 60601-1-2.

Regarding the static magnetic field, your MR system is continuously operated in the normal mode (means not more than 3.0 T). *For China only:* If the static magnetic field is higher than 2.0 T, the system is operated in the first level mode.

The limits of exposure of MR workers to static and time-varying magnetic fields may be regulated by local laws. Special precautions are necessary for pregnant MR workers, although no epidemiological evidence of any negative health effects currently exists (local laws may apply). It might be that the limits do not apply for pregnant MR workers. Additionally some countries may require a "member of the public" limit for pregnant MR workers, as the fetus is assigned as a member of the public.

6.5.4 Noise development

The exposure of MR workers to noise may be regulated by local laws. The operator must ensure that the sound level at the operating console is limited in compliance with local rules for the safety of the physician and the MR worker.

6.5.5 Laser

The laser of the laser light localizer is classified as Class 2M according to IEC 60825-1 (Class II according to US CDRH).

6.5.6 Video monitoring

Labeling obligation may be regulated by local laws. The operator is responsible to comply with these laws if a video system is installed.

6.5.7 Combination of devices

WARNING

Third-party devices connected to the MR system!

Several failures possible from electrical hazards, from RF burns to worsening of image quality

- ◆ Ensure that any combination with or modifications to the system comply with the requirements of IEC 60601-1 (chapter 16) and are accepted by Siemens.

Any application of physiological monitoring and sensing devices to the patient is done under the exclusive direction and responsibility of the system owner.

To perform a general test of the correct operation of peripheral equipment, a compatibility protocol is available. Manufacturers of peripheral equipment should use this Siemens compatibility protocol to test the functionality of its equipment.

The Physiological Measurement Unit (PMU) lets you control MR measurement sequences using a patient's physiological signals (ECG, respiration and pulse).

The PMU consists of the following components:

- PERU (Physiologic ECG and Respiratory Unit): ECG and respiratory sensor
- PPU (Peripheral Pulse Unit): Pulse sensor
- External trigger input

The physiological signals are acquired with receptors – ECG electrodes, respiratory cushion and pulse sensor - directly at the patient via the PERU (ECG, respiration) and PPU (pulse). The measured data can be viewed at the Dot display in the examination room and in the Physiological Display dialog window at the syngo Acquisition Workplace.

The Physiological Measurement Unit may be used only for controlling MR measurement sequences. The unit does not replace a patient monitoring system.

To prevent skin irritations, the PERU has to be located in the application cushion during the examination. The respiratory cushion is attached to the patient using the respiratory belt.

Both the PERU and the PPU are supplied with power via accumulators.

All other components of the PMU are supplied by system-internal voltage sources. The charging station is installed separately near the syngo Acquisition Workplace and is used for storing both units.

Accumulators should not be fully discharged before recharging them. If only one green LED flashes, you should charge the accumulator for the next patient. The maximum charging time is approx. 3 hours. After being fully charged, the operating hours for the units cover approximately 24 hours.

To charge a unit, it has to be placed firmly in the charging station. The PERU and PPU can be charged together or separately in the charging station.

6.5.8 Installation

⚠ WARNING Overload of circuit; water ingress; insufficient disconnection from the mains System malfunction; risk of electric shock ◆ Do not connect further components, which are not part of the original delivery volume, to the multiple socket of the computers. ◆ Do not connect further multiple sockets to the system. ◆ Ensure that the multiple socket of the computers is not placed on the floor. ◆ Do not use an isolating transformer.
⚠ WARNING Risk of electric shock ◆ To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

6.6 Emergency plans

Specific provisions must be taken for possible emergencies. This includes the generation of emergency plans (instructions on how to proceed/rescue scenarios) to prevent incorrect behavior under special circumstances.

6.6.1 Access to the examination room

Escape routes for the building must be established and well marked. Escape routes must not be obstructed.

- 1 Ensure that you have tools available to break the door open in case of an emergency.
- 2 Ensure that a window in the examination room can be used as an escape route in the case of an emergency.

⚠ WARNING RF door does not function as required! It is not possible to freely access or leave the examination room in case of an emergency ◆ Ensure that the RF door is checked and maintained regularly. ◆ Establish instructions on how to proceed in emergency situations. To avoid a situation where the patient is not accessible because the door cannot be opened (for example, the door handle becomes defective) establish alternative means/measures to enable the door to be opened quickly.

6.6.2 Emergency procedures

Instructions on how to proceed in emergency situations must be defined to ensure the safety of patients. For this purpose, MR-specific risks must be included (e.g. the presence of a magnetic field).

Among other things, the instructions must establish the fastest possible way for removing patients in emergency cases from the magnet's influence (if necessary, by shutting down the magnet). The instructions should also establish measures to ensure medical care as soon as possible.

As special precautionary measures, a program for medical supervision as well as a plan for using emergency equipment outside the magnet's influence must be in place for patients with a higher than normal risk factor, such as:

- Patients susceptible to cardiovascular collapse
- Patients who are likely to develop seizures
- Patients who are at an increased risk of heart attacks or other cardiac problems
- Patients with limited thermoregulation
- Claustrophobic patients
- Patients who are unconscious, anesthetized or confused or who are not able to communicate normally for other reasons
- Children (MR scanning is not established for safe imaging of fetuses and infants under two years of age. The responsible physician must evaluate the benefits of the MR examination compared to those of other imaging procedures.)

MR scanning is not established for safe imaging of fetuses and infants under two years of age. The responsible physician must evaluate the benefits of the MR examination compared to those of other imaging procedures.

6.6.3 Quench Emergency Plan

WARNING

Emergency situation with magnet quench!

Danger to life due to asphyxiation and hypothermia induced within a very short time

- ♦ Establish instructions on how to proceed in case gaseous helium escaping into the examination room.

The emergency plan should include the following information:

- Rescue scenarios that can be practiced with personnel
- Room-related conditions
- Rescue personnel (safety personnel, paramedics and firemen)
- Magnet Stop switch

Quench

During a quench, the super-conductivity of the magnet is lost. The energy of the magnetic field is converted into heat and the magnet field strength falls off. The liquid helium (coolant) boils off rapidly during this process and is released to the outside via the exhaust vent line. The escape of gaseous helium via the exhaust line is very noisy due to the high gas flow and you can hear a loud roaring noise for several minutes.

A quench may occur as follows:

- Start-up of the MR system (ramping up or filling the magnet)
- An accident (earthquake, fire, etc.)
- Spontaneously without any obvious external reason (highly unusual)

Defective vent line

If the vent line fails in part or fully, gaseous helium will enter the examination room. In this case the air conditioning unit will not be able to ensure sufficient air exchange and the following hazardous situations may arise:

- Poor visibility due to strong fog in the upper area of the room
- Rise in pressure in the examination room
- Hypothermia and risk of suffocation (e.g. in case of large leaks)

Due to such hazardous conditions as acute hypothermia and suffocation, rescue attempts must not be performed by a single person.

Persons not directly involved in the rescue should leave the examination room as well as adjacent rooms.

A filter (gas mask) without its own oxygen supply does not protect against suffocation through helium.

An emergency plan must be established to ensure correct conduct under such hazardous conditions.

6.6.4 Fire Fighting

In the event of fire, the fire has to be extinguished with methods appropriate to the surroundings. Respective firefighting equipment must be available. Fire precautions should be discussed with the local fire department and emergency procedures should be established. It is in the operator's responsibility to take the necessary initiatives.

Only fire extinguishers licensed for MR systems can be used.

Mandatory reporting in case of fire:

- Prior to initial start-up of the MR system, ensure that the fire department is informed about the MR system as well as structural on-site conditions.
- Inform the fire department about the contents of the measurement phantoms and the health risks caused by nickelous aerosol formations.

6.6.5 Employee qualifications/information

MR workers are individuals (e.g. operator, further personnel) who work within the controlled access area or MR environment. The system owner is responsible for ensuring that only trained and qualified MR workers and physicians are working on the MR system, so that they can perform all their tasks safely and efficiently, and in a way that minimizes their exposure to the electromagnetic field. In addition, the MR system may only be used as intended.

The system includes a key switch to prevent unauthorized switch on.

6.6.6 Pre-screening MR workers and patients

To lower the risks during exposure to the magnetic field, all patients as well as MR workers have to accomplish a pre-screening to avoid accidents and to establish safety measures. Therefore a pre-screening program shall be established by the operator, which helps the user to identify the patients and MR workers at risk. This especially applies to patients and MR workers who are at risk due to their professional activities, medical history and medical state as well as the influence of the MR equipment, such as:

1. Patients and MR workers:
 - with implants or with permanent make-up
 - with imbedded metal fragments from military activities
 - who are pregnant
2. Patients:
 - with typical contraindications
 - with a higher than normal likelihood of needing emergency medical treatment

6.6.7 Examination room and ambient conditions

6.6.7.1 Explosion protection

The MR system is not intended for operation in areas prone to explosive anesthetic gases.

6.6.7.2 Emergency switches

The voltage to the MR system can be turned off via a Power-Off (Emergency Shut-down) switch installed on-site. The switch can be used to switch off the electric power of the entire system immediately. The room installation must correspond to VDE 0100-710 and/or national laws.

A quench can be released with the Magnet Stop switch. This switch is installed on-site as well.

6.6.7.3 Air conditioning

An air conditioning system must be used to ensure the required environmental conditions

6.6.7.4 Magnetic fringe field and controlled access area

The fringe field can affect devices in the vicinity of the magnet. For this reason, the required safety distances must be observed. For details, please refer to the MR compatibility datasheet.

6.6.8 Signs and symbols

The system owner is responsible for properly identifying the accessible areas (e.g. regarding the electromagnetic field), the vicinity of the MR system, as well as adjacent areas by using the appropriate signs.

⚠ WARNING Missing hazard labels! Personal injury, property damage ◆ Attach the required warning and prohibition signs and observe national guidelines. ◆ Mark critical system areas with warning and prohibition symbols. ◆ Ensure that warning and prohibition signs are legible and clearly visible.

6.6.9 Maintenance/repair

6.6.9.1 Responsibility

As a supplier, Siemens will not be held responsible for the safety, reliability, and performance of the system in the following cases:

- Installations, additions, adjustments, modifications, and repairs to the MR system, or changes to the software that are not performed by Siemens Service.
- Assemblies are not replaced with original spare parts.
- The electrical wiring in the room does not meet the requirements of VDE regulation 0100-710 or applicable national regulations.

Siemens is not responsible for potential damage in the event non-authorized personnel refill the magnet with helium.

Modifications to the system or sub-systems during the service life require evaluation to the requirements of IEC 60601-1.

⚠ WARNING Unauthorized work on the magnet! Personal injury, property damage ◆ Only authorized personnel (Siemens Magnet Technology or Siemens) may perform work on the magnet. ◆ Do not open or remove safety valves and burst disks of the helium container. ◆ Do not change the standard configuration.
--

The monitors must only be opened by authorized Siemens Service personnel.

6.6.9.2 Maintenance at regular intervals

In the interest of the safety of patients, operating personnel, and third parties, it is strongly recommended that only authorized personnel perform the maintenance procedures provided by Siemens. System checks should be conducted more frequently if the system is operated under extreme conditions. For further information, see: System Owner Manual: Maintenance Plan.

Please inform Siemens Service if a maintenance contract does not exist.

In the course of system operation, technical and constructional changes may have been made to the MR system and its environment. It must be ensured that all components function satisfactorily and that the changes do not present hazardous conditions.

After system acceptance by the customer, a daily visual inspection should be performed for the MR system with respect to the following constructional changes:

- Changes in the environment at the output of the exhaust line (e.g. window installed retrospectively, inlet and outlet of air-conditioning units, new buildings or temporarily installed containers)
- Changes to the air conditioning unit or venting system (e.g. by adding air inlets and outlets in neighboring rooms)
- Installation of additional MR systems (e.g. prohibited use of the same exhaust line for several MR systems)
- Constructional changes inside and outside the examination room

The examination room must be checked for (newly introduced) magnetic parts.

The annual technical safety inspections are listed in this system owner manual and may only be performed by Siemens Service.

6.6.9.3 Repairs and modifications

All work, additions, and modifications to the MR system or to the installation site must be checked by Siemens in advance to ensure their compatibility with the MR system's functionality.

Modifications or additions to the product must comply with legal regulations.

The person performing the work must provide a certificate describing the nature and extent of work performed. This certificate must include information about changes to the nominal data or work area, along with the date, name of company, and signature.

Upon request, Siemens Service will provide technical documents for the MR system (e.g. circuit diagrams, spare parts lists, descriptions, calibration instructions). However, this does not constitute authorization for repairs.

6.6.9.4 Refilling helium

The magnet is filled with liquid helium as a coolant. Following installation, it is adjusted to the desired operating field strength.

During normal operation, the magnet does not lose helium. Under special conditions - power failure, malfunctions of the cold head and maintenance activities - liquid helium must be refilled by Siemens Service or specially qualified personnel.

When filling the magnet with helium, perform the necessary tasks carefully and accurately, observing all regulations. Wear protective clothing to prevent frostbite.

Liquid helium presents the following properties that, among other things, may result in hazardous conditions when not handled professionally:

- Extremely cold: causes frostbites when it comes in touch with skin
- Oxygen in ambient air is displaced during boil-off: risk of asphyxiation

If the helium fill level is too low, the alarm box or the syngo Acquisition Workplace will signal this accordingly. In case of alarm, notify Siemens Service and/or ensure refilling only through trained and experienced personnel.

 **WARNING**

Escaping gaseous helium during refill activities!

Injury to persons, danger of suffocation, frostbite

- ◆ Ensure that the rooms are ventilated via an air conditioning system. This includes refilling with helium as well.
- ◆ Ensure that escape routes have been determined, are identified as such and are not obstructed.
- ◆ Ensure that the magnet is only filled by Siemens Service.

It is prohibited to store flammable material in the vicinity of containers filled with coolant.
Use non-magnetic coolant containers for the helium.

7. Environment Protection Requirements

Siemens Healthcare GmbH has a certified environmental management system, as meeting the requirements of ISO 14001:2004 for the research, development, manufacture, assembly, sale, supply, service, administration, training, serviceability, documentation and refurbishment of diagnostic imaging systems.

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9. The list of International Standards / Directives / Documents applicable to the Medical Device (ISO, CE, DoC, IEC)


IEC 60601-1:2012 EN 60601-1:2006 + Cor.:2012 + A1:2013	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
IEC 60601-2-33:2010/A1:2013 /A2:2015 EN 60601-2-33:2010/A2:2015	Medical electrical equipment - Part 2-33: Particular requirements for the basic safety and essential performance of magnetic resonance equipment for medical diagnosis
IEC 60601-1-2:2014 EN 60601-1-2:2015	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances - Requirements and tests
IEC 60601-1-6:2010 EN 60601-1-6:2010	Medical electrical equipment – Part 1-6: General requirements for safety – Collateral standard: Usability
IEC 60601-1-6:2010+A1:2013 ¹ EN 60601-1-6:2010+A1:2015 ¹	Medical electrical equipment – Part 1-6: General requirements for safety – Collateral standard: Usability
IEC 60825-1:2014 EN 60825-1:2014	Safety of laser products - Part 1: Equipment classification, requirements and user's guide
IEC 62304:2006+A1:2015 EN 62304:2006+Cor.:2008+ A1:2015	Medical device software - Software lifecycle processes
IEC 62366:2007 EN 62366:2008	Medical devices - Application of usability engineering to medical devices
IEC 62366-1:2015 ¹ EN 62366-1:2015 ¹	Medical devices - Part 1: Application of usability engineering to medical devices
EN 1041:2008/A1:2013	Information supplied by the manufacturer of medical devices
EN ISO 15223-1:2016	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
ISO 10993-1:2009 EN ISO 10993-1:2009+AC2010	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
ISO 14971:2007 EN ISO 14971:2012	Medical devices - Application of risk management to medical devices



¹ For MAGNETOM Aera with syngo MR E11C-AP01 and MAGNETOM Systems with syngo MR E11E only

10. List of equipment and measurement tools necessary for quality control of the medical device

10.1 Measurement Phantoms




The following accessories for quality assurance are available at the customer.


	Plastic bottle, 15ml ■ Solutes per 1000 g H₂O dist.: 3g MnCl₂ * 4 H₂O + 5g NaCl
	Plastic bottle, 1900ml ■ Solutes per 1000 g H₂O dist.: 3.75g NiSO₄* 6 H₂O + 5g NaCl
	Plastic bottle, 5300ml ■ Solutes per 1000 g H₂O dist.: 3.75g NiSO₄* 6 H₂O + 5g NaCl
	Spherical phantom D165 ■ Diameter: 165 mm ■ Contents: 2570 ml ■ Solutes per 1000 g H₂O dist.: 1.25g NiSO₄ * 6 H₂O
	Spherical phantom D240 ■ Diameter: 240 mm ■ Contents: 7300 ml ■ Solutes: MARCOL-Oil + 0.011g MACROLEX blue

	Hand/Wrist 16 phantom ■ Contents: 1100 ml ■ Solutes per 1000 g H₂O dist.: 3.75g NiSO₄* 6 H₂O + 5g NaCl
	Foot/Ankle 16 phantom ■ Contents: 3300 ml ■ Solutes per 1000 g H₂O dist.: 3.75g NiSO₄* 6 H₂O + 5g NaCl
	Special phantom for Endorectal A powder is included with the special phantom for the Endorectal coil. It is used to prepare the phantom fluid. For additional information, refer to the Operator Manual for the special phantom.

10.2 Phantom Holder

The following phantom holders need to be used for the coil's quality assurance.

 Phantom holder 17	 Phantom holder 18	 Phantom holder 23
 Phantom holder 24	 Phantom holder 25	 Phantom holder 26

11. Transportation requirements.

11.1. Transportation requirements.

- The magnet can be shipped and transported worldwide with or without liquid helium. Maximum possible helium level during shipping: 95%.
- The magnet is classified as "Dangerous goods" as defined by the ADR, IMDG and IATA codes and their requirements. Liquid helium is classified under UN1963 (Helium, refrigerated liquid).
- Transportation of the system containing liquid helium must be in accordance with the appropriate national and international codes.

- Using an air-ride truck with or without a shipping pallet the following maximum peak acceleration measured at the magnet mounting points during transportation shall not exceed:
 - 5.0g vertical shock load applied to the cryostat mounting points.
 - 2.0g horizontal shock load applied to the cryostat mounting points in any direction.
 - Shock is defined as a dynamic force with duration of no more than 0.5 seconds.
- Road transport is allowed up to 15.000 miles.
- Road transport may be with or without an approved shipping pallet. Without NO shipping pallet, the magnet must be supported using transport beams.
- Road vehicles must be fitted with air ride suspension, whether or not a shipping pallet is used.
- Always transport a magnet using vehicles fitted with air ride suspension.
- Road vehicles must be suitable for the intended usage, be properly maintained at all times, and comply with all local laws and regulations.
- Relevant requirements of the ADR code shall be satisfied.
- The vehicle and operator must be approved by Siemens.
- This specification does not constitute acceptance for use as part of a mobile system. Mobile magnets require a specific trailer design for continuous road transportation.
- Air transport is allowed up to 30.000 miles.
- If air transport is used for any section of the route, an approved shipping pallet must be used for the entire transportation process.
- Relevant requirements of the IATA code must be satisfied.
- Sea transport is allowed up to 30.000 miles.
- If sea shipment is used for any section of the route, an approved (magnet system related) shipping pallet must be used for the entire transportation process.
- Relevant requirements of the IMDG code must be satisfied.
- Rail Transportation is not allowed!
- Transportation through the channel tunnel (UK to France) is permissible if the system is mounted on an approved air ride truck. The air ride suspension may be inactive.
- Do not use a fork lift for lifting and transporting the magnet!
- Lifting the magnet must always be conducted by personnel who are trained in handling components and medical equipment of this size and weight.
- Lifting during transportation including the final movement onto the installed site must be performed using a suitable sized crane capable of lifting the full magnet weight.

- Suitable straps or chains capable of lifting the full magnet weight must be used. The straps or chains must be secured to the magnet lifting points using the four shackles supplied and fitted to the magnet.
- Use four suitable rated slings or chains (minimum of 3 m in length) for lifting the full system weight.
- Check the magnet components are not fouled by the lifting straps or chains.
- Do not stand on the Magnet. Damage to components may occur.
- All transport activities including delivery of the magnet and the system components to the MR suite may be performed by a licensed carrier only.

11.2. Transportation specifications – Magnet

Transport specifications for the Magnet	
Storage temperature	-20 °C to +60 °C
Relative humidity	10% to 90% non-condensing (must be kept dry if unpacked)
Atmospheric pressure	20 kPa to 106 kPa

11.3. Transportation specifications – System Components

Transport specifications for the system components	
Ambient temperature	-25 °C to +60 °C
Relative humidity	10% to 95%
Atmospheric pressure	70 kPa to 105 kPa

12. Packaging

The dimensions of the packaging are Length 350 cm, Width 215 cm, Height 235 cm.

The magnet including cover has to be handled with care and has to be fixed in the custom-made transport box, all packaging labels and safety labels are located on the packaging. Please refer to picture below.

The electronic cabinet (ACC, GPA XXL, cables technic-room -> operator-room) is packed in a wooden transport box, all packaging labels and safety labels are located on the packaging. Please refer to picture below.

The computer, monitor and accessories are packed in a cardboard transport box, all packaging labels and safety labels are located on the packaging. Please refer to picture below.

The Filter-Panel (including cables filter-panel -> ACC) is packed in a wooden transport box, all packaging labels and safety labels are located on the packaging. Please refer to picture below.

The System Separator (incl. He-Compressor) is packed in a wooden transport box, all packaging labels and safety labels are located on the packaging. Please refer to picture below.

The Helium Compressor is packed in a cardboard transport box, all packaging labels and safety labels are located on the packaging. Please refer to picture below.

The external hose kit is packed in a cardboard transport box, all packaging labels and safety labels are located on the packaging. Please refer to picture below.

Test phantoms, phantoms-holder and accessories for patient positioning are packed in a cardboard transport box, all packaging labels and safety labels are located on the packaging. Please refer to picture below.

Look-Covers are packed in a wooden transport box, all packaging labels and safety labels are located on the packaging. Please refer to picture below.

The cable set (Gradient cables, Ramp cables, Accessories) is packed in a cardboard transport box, all packaging labels and safety labels are located on the packaging. Please refer to picture below.

The surface coils and documentation are packed in a cardboard transport box, all packaging labels and safety labels are located on the packaging. Please refer to picture below.

Patient table is packed in a wooden transport box, all packaging labels and safety labels are located on the packaging. Please refer to picture below.

13. Labeling

Packaging Label for MAGNETOM Aera until mid of July 2016 was as following:

Packaging Label for MAGNETOM Aera starting mid of July 2016 is as following:

14. Storage and Service Life

Siemens considers the subject devices exempt from having a shelf life due to the fact that the devices are manufactured upon order and is therefore not susceptible to degradation that would lead to functional failure.

Once an order is placed and a machine is manufactured, it is transported, installed and QC and tune up tests are completed on site to ensure proper function of the MR scanner. The MR scanners do not spend time on the "shelf" that would pose a risk to the performance of the device or its components.

The expected service life of 10 years of operation under the environmental conditions indicated in the instructions for use as well as with all preventive and corrective maintenance performed by qualified technical personnel authorized by Siemens Healthcare GmbH.

If the systems are used beyond this time, additional maintenance and service may be required in order to maintain system performance and reliability.

14.1 Storage specifications – Magnet

Storage specifications for the Magnet	
Storage temperature	-20 °C to +60 °C
Relative humidity	10% to 90% non-condensing (must be kept dry if unpacked)
Atmospheric pressure	20 kPa to 106 kPa

14.2 Storage specifications – System Components

Storage specifications for the system components	
Ambient temperatur	-25 °C to +60 °C
Relative humidity	10% to 95%
Atmospheric pressure	70 kPa to 105 kPa

15. Procedure of recycling and disposal (Medical Device recycling or disposal data)

On disposing of the complete system or parts thereof, currently valid environmental legislations must be observed. Examples of environmentally relevant components are:

- Accumulators and batteries;
- Transformers;
- Capacitors.

All medical waste shall be collected, used, decontaminated, placed, stored, transported, listed and disposed of in accordance with the laws on public health and epidemiological welfare.

System components hazardous to persons or the environment must be disposed of with care and in compliance with legally binding ordinances.

For details the local customer service representative shall be contacted.

16. Usage conditions

16.1. Installation requirements

To ensure correct installation, complete that all preliminary steps specified. For detailed information on preparing the room, transporting the system, and other installation requirements, consult the project manager. The on-site project manager is responsible for all room and system-related preparation.

The following prerequisites must be fulfilled:

- Prepare and clean the room.
- Make sure the power and water connections are available.
- The internal resistance of the on-site power supply unit must correspond to the requirements.
- The transport route for the system must be defined and free of obstructions (regarding dimensions, access, load).
- Ensure that the magnet room is properly equipped.
- Ensure that the openings of the RF shield are sufficiently large.
- Ensure that there is sufficient space to install the magnet.
- Take account of narrow passages due to suspended ceilings, carriers, and columns.
- Determine the final location of the magnet in the magnet room.
- Ensure that the primary cold water supply is ready for operation.
- Block off the unloading zone to ensure that passersby are not harmed.
- Do not block escape routes or essential access routes to the hospital.
- Equip the rooms properly. The RF-shield, including the opening for inserting the magnet, must be installed. The magnet room must also be equipped with sufficient lighting.
- The floor must be clear so that the magnet and system components can be transported.
- The four corners of the magnet must be drawn on the floor of the magnet room to determine the magnet's location. Refer to the project plan.
- The floor from the unloading zone to the final location must be covered with metal plates (at least 5 mm thick).
- The tools required for installation must be available on-site (steel rollers, hydraulic, lifts, etc.).

16.2 Information about maintenance (general and on the territory of Russia)

Routine checks and periodic maintenance are necessary to ensure safe and satisfactory operation of the system. Proof of these activities is required by authorities in certain countries.

Routine checks include:

Daily, weekly, and monthly inspections as well as legally required checks as described in the system "Function and Safety Checks" chapter of the operating instructions. Normally, the system operator entrusts the clinical operating personnel with the task of performing these routine checks.

Periodic maintenance includes:

- Safety check (including safety-related tests)
- Preventive maintenance
- Quality and function tests
- Replacement of safety-relevant parts subject to wear.

This work has to be performed by qualified and authorized service engineers only. In this context, qualified means that the engineers have been trained accordingly or have acquired practical experience through routine service activities. Authorized means that the engineers have been authorized permission by the operator of the system to perform maintenance work.

Upon first start-up of the system, Siemens ask the system operator to designate a staff member responsible for ensuring that routine checks, preventive inspection, and maintenance work are performed. This staff member is responsible for archiving all certificates in the "System Owner Manual" binder.

In addition to our repair service, Siemens also offers a complete range of services for the preventive inspection and maintenance of the system. These services can be called on as required or agreed upon in a flexibly drafted maintenance contract.

Service Contacts are:

Siemens Healthcare GmbH – Customer Service (CS)
Address: Henkestr. 127, 91052 Erlangen, Germany
Telephone: +49 9131 84-0
Web: www.siemens-healthineers.com

Siemens Healthcare LLC
Address: Dubininskaya Str. 96, 115093 Moscow, Russia
Tel.: +7 495 737 12 52
Web: www.siemens-healthineers.com/ru

16.3 Information about calibration of the MD / requirements to calibration (if applicable)

There is no calibration required for the MR systems. Nevertheless, there are different checks and tests which have to be done as described in the system owner manual.

17. Warranty

The system owner is responsible for the software licenses including software warranty conditions.

The specific warranty conditions are usually defined in the contractual agreement between Siemens Healthcare GmbH and the system owner.

18. Protective precautions for use

The MAGNETOM MR system must be operated according to the intended use and only by qualified persons with the necessary knowledge in accordance with country-specific regulations, e.g. physicians, trained radiological technicians or technologists, subsequent to the necessary user training.

This user training must include basics in MR technology as well as safe handling of MR systems. The user must be familiar with potential hazard and safety guidelines the same way the user is familiar with emergency and rescue scenarios. In addition, the user has to have read and understood the contents of the operator manuals.

The system owner manual and operator manuals contain the most important safety aspects for which the owner/operator of the MR system are responsible. These include legal requirements, emergency plans, employee information and qualifications, as well as requirements that must be met in the examination room.

Operating personnel as well as personnel who are not regulars in the examination room (e.g. cleaning personnel, rescue personnel) must be informed about the special conditions in the examination room. Magnetizable devices (e.g. floor polishers, vacuum cleaners, wheelchairs, metal gurneys) must not be used in the examination room.

19. Contact for reclamation

In case of any questions or reclamations, the following shall be contacted:

Germany

Siemens Healthcare GmbH - Magnetic Resonance (MR)
Address: Henkestr. 127, 91052 Erlangen, Germany
Telephone: +49 9131 84-0
Web: www.siemens-healthineers.com

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Herewith, I certify the foregoing signature,
acknowledged before me, of

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Mrs. Friederike Hertle is personally
known.

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CE 0123

Insertion Plate Breast Coil 1.5T		Medical equipment for MRI imaging	
Manufacturer: NORAS MRI products GmbH Leibnizstr. 4 97204 Höchberg Germany		Distributed by: Siemens Healthineers GmbH	
Series Ref.: 112903		Rev. 06	
		Rev. 06	
www.noras.de/insert-plate			

Insertion Plate Breast Coil 1.5T / 63,6 MHz

Manufacturer: Noras MRI products GmbH Leibnizstrasse 4 97204 Höchberg Germany	Distributor: Siemens AG Wittelsbacher Platz 2 80333 München Germany
--	--

CE 0123

Patient Pad Breast Coil 1.5T		Medical equipment for MRI imaging	
Manufacturer: NORAS MRI products GmbH Leibnizstr. 4 97204 Höchberg Germany		Distributed by: Siemens Healthineers GmbH	
Series Ref.: 112902		Rev. 07	
		Rev. 07	
www.noras.de/patient-pad			

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


Rev (2P) 06

Distributed by:
Siemens AG
Wittelsbacherplatz 2
DE-80333 Muenchen
Germany

100% Subcontracted
Manufactured in
the USA
Invivo Corporation
Gainesville, FL 32608-7991
U.S.A. +1 (352) 336-0010

	 16 CH BREAST SEN MR COIL KIT 1.5T REF  (1P) 10 60 66 09 SN  (S) serial Rev  (2P) 03 Distributed by: Siemens Healthcare GmbH Henkestr. 127 91052 Erlangen Germany CE 0123 EC REP GBM Authorized Representative Ltd The White House 2 Meadrow Goswilling GU7 3HN United Kingdom +44 (0)7710 039 721 www.mba-gbm.com
	 16 CH BREAST SEN MR COIL KIT 1.5T Distributed by: Siemens AG Wittelsbacherplatz 2 DE-80333 Muenchen Germany CE 0123 GBM Authorized Representative Ltd The White House 2 Meadrow Goswilling GU7 3HN United Kingdom +44 (0)7710 039 721 REF  (1P) 10 00 05 09 SN  (S) serial Rev  (2P) 03
Катушки для визуализации и биопсии молочных желез (Breast BI 7 Coil)	Breast BI 7 MR Coil 1,5T GTIN Origin Rev Serial Model (81)8425129920259 (42)276 (20)81 (71)XXXX (248)18676741 REF 118879 SN XXXX 2019-07-08 Rev. 01 Distributed by: Siemens Healthcare GmbH CE 0123
Катушки эндоректальные для визуализации и спектроскопии простаты	 PROSTATE 2 SEN MR COIL KIT 1.5T REF  (1P) 10 60 66 24 SN  (S) serial Rev  (2P) 002 Distributed by: Siemens Healthcare GmbH Henkestr. 127 91052 Erlangen Germany CE 0123 EC REP GBM Authorized Representative Ltd The White House 2 Meadrow Goswilling GU7 3HN United Kingdom +44 (0)7710 039 721 www.mba-gbm.com

	   PROSTATE 2 SEN MR COIL KIT 1.5T REF:  (1P) 10 60 66 24 SN:  (S) XXXXXX Rev:  (2P) 002 Distributed by: Siemens AG Wittelsbacherplatz 2 DE-80333 Muenchen Germany For Rx Use Only 
Катушки для небольших участков тела	4 - Channel Special - Purpose Coil  (01) 04056869014968 GTIN (240) 10499077 MODEL (422) 156 ORIGIN CN (21) 1625 SERIAL IVK Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center, Gaoxin C. Ave., 2nd, Hi-Tech Industrial Park, 518057 Shenzhen, China Made in China  (1P) MODEL No. 10499077 (S) SERIAL/LOT No. 1017 MADE IN CHINA PS 02 03 04 05 06 07 08 09 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24 25 26 27 28 MR2 AAN330071 
Катушки гибкие большие, дополнительные (Flex Large 4 Coil)	SIEMENS Flex Large 4 MR coil 1.5T  (01) 04056869014784 GTIN (240) 10185552 MODEL (422) 276 ORIGIN DE (21) 30936 SERIAL IVK Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany Made in Germany ES: 02 03 04 05 06 07 08 09 10 11 12 13 14 15 16 17 18 19 SIEMENS (1P) Model No. 10185552 (S) Serial No. 4182 IVK Siemens AG, Wittelsbacherplatz 2, DE-80333 Muenchen Made in Germany 

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IEC 60601-1:2012 EN 60601-1:2006 + Cor.:2012 + A1:2013	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
IEC 60601-2-33:2010/A1:2013 /A2:2015 EN 60601-2-33:2010/A2:2015	Medical electrical equipment - Part 2-33: Particular requirements for the basic safety and essential performance of magnetic resonance equipment for medical diagnosis
IEC 60601-1-2:2014 EN 60601-1-2:2015	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances - Requirements and tests
IEC 60601-1-6:2010 EN 60601-1-6:2010	Medical electrical equipment – Part 1-6: General requirements for safety – Collateral standard: Usability
IEC 60601-1-6:2010+A1:2013 ¹ EN 60601-1-6:2010+A1:2015 ¹	Medical electrical equipment – Part 1-6: General requirements for safety – Collateral standard: Usability
IEC 60825-1:2014 EN 60825-1:2014	Safety of laser products - Part 1: Equipment classification, requirements and user's guide
IEC 62304:2006+A1:2015 EN 62304:2006+Cor.:2008+ A1:2015	Medical device software - Software lifecycle processes
IEC 62366:2007 EN 62366:2008	Medical devices - Application of usability engineering to medical devices
IEC 62366-1:2015 ¹ EN 62366-1:2015 ¹	Medical devices - Part 1: Application of usability engineering to medical devices
EN 1041:2008/A1:2013	Information supplied by the manufacturer of medical devices
EN ISO 15223-1:2016	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
ISO 10993-1:2009 EN ISO 10993-1:2009+AC2010	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
ISO 14971:2007	Medical devices - Application of risk management to medical devices

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