Medtronic

18000 Devonshire Street
Northridge, CA 91325 USA
www.medtronic.com
Tel 818-576-5211

ATTESTATION STATEMENT

As a responsible member of Medtronic MiniMed, I hereby attest that I believe the Instructions for Use attached herein, to be a true and accurate copy.

Katherine Rojas
Regulatory Affairs Specialist

Kath Rojas July 26, 2021
Signature

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document


State of California
County of Ventura

On July 26th 2021 before me, Kusum Manwani Notary Public
(date) (insert name and title of the officer)

personally appeared Katherine Rojas, who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the same in his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the person(s), or the entity upon behalf of which the person(s) acted, executed the instrument. I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

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
Medtronic

© 2019 Medtronic. Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic.

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Contacts:

Africa:

Medtronic South Africa and Southern Africa

Office Reception Tel:
+27(0) 11 260 9300

Diabetes: 24/7 Helpline:
0800 633 7867

Sub-Sahara 24/7 Helpline:
+27(0) 11 260 9490

Albania:

Net Electronics Albania
Tel.: +355 697070121

Argentina:

Corpomedica S.A.
Tel.: +(11) 4 814 1333
Medtronic Directo 24/7: +0800 333 0752

Armenia:

Exiol LLC
Tel.: +374 98 92 00 11
or +374 94 38 38 52

Australia:

Medtronic Australasia Pty. Ltd.
Tel.: 1800 668 670

Bangladesh:

Sonargaon Healthcare Pvt Ltd.
Mobile: (+91)-9903995417
or (+880)-1714217131

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USA:

Medtronic Diabetes Global
Headquarters
24-Hour Technical Support:
+1-800-646-4633
To order supplies: +1-800-843-6687

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Introduction

The Quick-serter (MMT-305QS) is used as an aid to insert compatible Medtronic MiniMed infusion sets.

Indications for use

The serter is indicated for single-patient use. It is not intended for multiple patients. Refer to the system user guide for additional information, such as the intended age group.

Contraindications

No contraindications are associated with serter use. For contraindications related to continuous glucose monitoring, see the system user guide.

Clinical benefits

Not applicable. The serter is an accessory that helps to insert compatible Medtronic infusion sets.

User safety

Warnings

Read this entire user guide prior to inserting the infusion set. Failure to follow directions may result in improper insertion, pain, or injury.

Improper insertion can result in non-delivery of insulin, which can lead to hyperglycemia.

Always inspect the infusion set packaging for damage before use. Infusion sets are sterile and non-pyrogenic unless the package has been opened or damaged. If the infusion set packaging is open or damaged, dispose of the infusion set according to local regulations. Use of a non-sterile infusion set can result in infection.

Failure to secure the infusion set in the serter prior to insertion, with the needle pointing straight up, may result in improper insertion, pain, or minor injury.

Never point the loaded serter toward any body part where insertion is not desired.

Do not use the infusion set if the introducer needle is bent. Always inspect the introducer needle after the needle guard has been removed. If the introducer needle is bent, discard the infusion set.

For questions or concerns related to serter use, contact a local Medtronic support representative for assistance.

For medical questions or concerns, or to report a serious incident, contact a healthcare professional.

Precautions

The first use of the Quick-serter should be done in the presence of a healthcare professional. Improper insertion may lead to infection, diabetic ketoacidosis, or hyperglycemia.

Risks and side effects

Risks and side effects include infection, diabetic ketoacidosis (vomiting, dyspnea, nausea), and hyperglycemia (polydipsia, frequent urination).

Hazardous substances

None.

Allergens

None known.

Maintenance

Cleaning

To clean the serter:

1. Rinse the serter under running tap water.
2. Submerge the serter in a solution of mild soap and warm water.
3. Use a soft-bristled brush to clean the serter until it is free of visible dirt.
4. Rinse the serter under running tap water until all visible soap is gone. Make sure to fully rinse all hard-to-reach areas.
5. Shake off any excess water. Place the serter upright on a clean cloth to air dry.

Storage

Store the serter in the released position at room temperature.


Disposal

Dispose of the serter according to local regulations.

Warranty

For product warranty information, please contact a local Medtronic support representative.

Components

- (A) Handle
- (B) Tubing slot
- (C) Release button
- (D) Side buttons

Instructions

Unfold the front cover of this booklet to view images that are numbered with the steps below. Refer to the infusion set user guide before using the Quick-serter. The first use of the Quick-serter should be done in the presence of a healthcare professional.

Preparation

1. Wash hands with soap and water.
2. Clean the insertion site. Use an alcohol wipe or other disinfectant as directed by a healthcare professional. The area must be dry before the infusion set is inserted.

Using the Quick-serter

3. Place the Quick-serter onto a sturdy flat surface with the handle facing downward. Place the infusion set into the Quick-serter. If needed, place infusion set tubing in the tubing slot.

4. Place the infusion set securely within the Quick-serter. To do this, place fingers on the adhesive liner and push the infusion set down.
5. Remove both pieces of the adhesive liner.
6. Pick up the Quick-serter and pull the handle down until it clicks.
7. Twist the needle guard to loosen it. Carefully remove the needle guard. Do not touch or bend the needle.
8. Place the Quick-serter on the desired infusion site. Then, press the two side buttons at the same time to insert the infusion set.
9. Press the release button on top of the Quick-serter to release the infusion set.
10. Hold the infusion set tubing with one hand. Then, remove the Quick-serter by lifting it away from the body with the other hand.
11. Pull the needle housing away from the infusion set with one hand. Gently press the tape against the skin with the other hand.
12. Close the needle housing. Then, discard the shielded needle in a sharps container.
13. Smooth out the tape securely against the skin.

Technical specifications

Device specifications	
Storage Conditions	Temperature range of 20 °C (68 °F) to 30 °C (86 °F)
	Relative humidity range of 15% to 93%

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ATTESTATION STATEMENT

As a responsible member of Medtronic MiniMed, I hereby attest that I believe the following Appendix to the IFU attached herein, to be a true and accurate copy.

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
On March 11th 2022 before me, Kusum Manwani Notary Public
(date) (insert name and title of the officer)

personally appeared Katherine Rojas, who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the same in his/her/their authorized capacity(ies); and that by his/her/their signature(s) on the instrument the person(s), or the entity upon behalf of which the person(s) acted, executed the instrument. I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

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Sil-Serter (MMT-385)

The Sil-serter (MMT-385) is manufactured, packaged, and labeled by SMC Ltd per Medtronic specifications. It is a non-sterile, standalone device packaged in a cardboard box of pre-defined dimensions along with folded copies of IFU. A product identifying label is affixed to the side of the box. The product is distributed to the patients in the primary packaging and does not have any additional secondary packaging.

Contents of delivery:

- Sil-serter, model MMT-385 (1 unit)
- Instructions for Use

Product in the amount of 168 pieces are placed in a cardboard transport box with "bubble" and sealed with packaging tape.

Devices will be used in the following operating conditions:

Quick-Serter:

Operating Conditions

Temperature, °C	+5 to +40
Humidity, %	15 – 93, non-condensing
Pressure, kPa	70 - 106

Sil-Serter:

Transport, Use and Storage Conditions

Temperature, °C	-20 to +50
Humidity, %	Up to 90
Pressure, kPa	N/A

A sterile use space is not required to operate the serters. However, preparation steps such as cleaning the insertion site are required. Refer to the instructions for use for more information.

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Read all instructions carefully before using. For information about pump therapy, potential errors, and possible risks and their causes, consult your pump user guide before connecting the infusion set. Refer to your infusion set user guide before using the. If you are using for the first time, you should perform your first set-up in the presence of your healthcare professional.

Insertion device has been designed to minimize the risk of accidental needle sticks and inadvertent actions such as "launching" the infusion sets; however, users are warned that improper use of the device may lead to insertion pain or user injury.

The device has a double button ejection to prevent accidental discharge; both buttons must be pressed to release the spring/carrier from the locked position. The device also has a release button for positive disengagement of the infusion set from the insertion device.

Always inspect the packaging of the infusion set for damage. Do not use the infusion set if there is damage.

Never point the loaded toward any body part where insertion is not desired. This is to avoid an accidental needle stick.

Carefully remove the needle guard before inserting the infusion set into your skin. The needle should be straight and point outward. Do not use the infusion set if the needle is bent.

You should perform your first set-up of in the presence of your healthcare professional.

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Medical devices will be stored in a warehouse environment after the product is manufactured, and before delivery to the customer. The transportation must be done in the covered vehicles of all types in accordance with the rules in force in used type of transport.

Quick-Serter:

Storage Conditions (home use)

Temperature, °C	Room Temperature (+20 °C to 30 °C)
Humidity, %	15 to 93
Pressure, kPa	N/A

Transport and Storage Conditions (warehouse)

Temperature, °C	-25 to +70
Humidity, %	15 to 93

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