Medtronic

Medtronic
Medtronic Xomed, Inc.
6743 Southpoint Drive North
Jacksonville, Florida 32216
USA

Custodian:

Alexandra Oliver

Alexandra Oliver
Assoc. Regulatory Affairs Specialist
Medtronic Xomed Inc.

12/12/2017
Date (MM/DD/YYYY)

Instructions For Use

XPS Nexus System & Accessories

ATTESTATION

STATE OF FLORIDA
DUVAL COUNTY

On this 12th day of December 2017, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following, Instructions For Use, issued by Medtronic Xomed, Inc., presented to me by the document's custodian, Alexandra Oliver, and, to the best of my knowledge, that the photocopied document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

Notary:

Inez R. Anderson
Official Notary Signature

Subscribed and sworn to before me
This 12th day of December, 2017.

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Medtronic

Medtronic
Medtronic Xomed, Inc.
6743 Southpoint Drive North
Jacksonville, Florida 32216
USA

Custodian:

Alexandra Oliver

Alexandra Oliver
Assoc. Regulatory Affairs Specialist
Medtronic Xomed Inc.

12/12/2017
Date (MM/DD/YYYY)

Instructions For Use

XPS Nexus System & Accessories

ATTESTATION

STATE OF FLORIDA
DUVAL COUNTY

On this 12th day of December 2017, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following, Instructions For Use, issued by Medtronic Xomed, Inc., presented to me by the document's custodian, Alexandra Oliver, and, to the best of my knowledge, that the photocopied document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

Notary:

Inez R. Anderson
Official Notary Signature


Subscribed and sworn to before me
This 12th day of December, 2017

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Table of contents

GLOSSARY	4
INDICATIONS FOR USE	5
DEVICE DESCRIPTION	5
CONTRAINDICATIONS	5
WARNINGS	5
PRECAUTIONS	6
WHEN THE SYSTEM ARRIVES	7
ATTACHMENT ASSEMBLY	7
IRRIGATION ASSEMBLY	7
TOOL OPERATION	7
SYSTEM LIGHTS AND FIGURES	7
System lights	7
System figures	8
SETUP THE XPS NEXUS	13
CHANGE THE DEFAULT SETTINGS ON THE XPS NEXUS	13
CONSOLE REUSE AND CLEANING	14
STERILITY	14
SERVICE AND MAINTENANCE	14
INDIGO HIGH-SPEED OTOLOGIC DRILL	14
Device description	14
Indigo drill straight attachment assembly	14
Indigo drill angled attachment assembly	15
Indigo drill irrigation and tubing clip assembly	16
Indigo drill cleaning and sterilization instructions	16
Indigo drill technical specifications	16
STRAIGHTSHOT M4 AND STRAIGHTSHOT MAGNUM II MICRODEBRIDERS	16
Device description	16
StraightShot M4 blade or bur assembly	17
StraightShot M4 suction and irrigation tube assembly	17
StraightShot Magnum II blade and bur assembly	18
Microdebrider blade control	18
Cleaning and sterilization instructions	18
StraightShot M4 technical specifications	18
StraightShot MAGNUM II technical specifications	18

TROUBLESHOOTING..... 19
 XPS Nexus and foot control 19
 MEDTRONIC XOMED blades or burs 19
SYSTEM REQUIREMENTS AND SPECIFICATIONS..... 20
 Console specifications 20
GUIDANCE AND MANUFACTURER’S DECLARATION – ELECTROMAGNETIC IMMUNITY 21
 Part I 21
 Part II 22
RETURNS 23
CUSTOMER SERVICE INFORMATION 23
LIMITED WARRANTY 23
REPROCESSING INSTRUCTIONS..... 24
SYMBOLS 44

Glossary

The following words and acronyms may be used in this guide.

CC/min	Flow rate
FCU	Foot Control Unit
FWD	Forward - Rotation is clockwise
I.V.	Intravenous
OSC	Oscillate
REV	Reverse - Rotation is counter-clockwise
RPM	Rotations per minute
XPS	Xomed Power System

EN**Indications for use**

The XPS Nexus System is indicated for the incision/cutting, removal and drilling of soft and hard tissue and bone in head and neck/ENT, oral/maxillofacial and plastic/reconstructive/aesthetic surgical procedures.

Device description

The XPS Nexus System consists of a power control console, foot control, connection cables, and assorted handpieces to drive various burs, blades, drills, rasps, and cannulae.

The XPS Nexus console provides power to the attached handpiece. In addition to the StraightShot M4 Microdebrider the console also powers the StraightShot Magnum II and Indigo handpieces. It includes an integrated irrigation pump for irrigation of blades and burs.

The StraightShot M4 Microdebrider features a finger wheel that enables the user to rotate the cutting tip of specially designed straight and curved rotatable blades. It is also compatible with non-rotating blades and burs. It includes a locking mechanism to prevent inadvertent rotation of non-rotating blades. The lightweight design combines high power and precision performance and can operate in forward mode of 12000 rpm for bur operation and 5000 rpm in oscillate for blade operation. The StraightShot M4 also provides irrigation tubing management grooves to keep tubing out of the surgeon's way during a procedure.

The StraightShot Magnum II handpiece's lightweight design combines high power and precision performance and can operate in forward mode of 12000 rpm for bur operation and 5,000 rpm oscillate for blade operation.

The Indigo handpiece is a small, compact high-speed, high-torque, reversible electric drill that can be used to dissect bone and biomaterial at variable speeds up to 60,000 rpm.

This device is intended for use by physicians trained in the procedures described, in any healthcare facility environment.

Medtronic attachments are biocompatible per recommendations set forth in ISO 10993-1.

Part Description	CFN Number
XPS Nexus Console	1894510
XPS Nexus Foot Control	1894520
Indigo high-speed otologic drill	1845000
StraightShot Magnum II Microdebrider	1897200
StraightShot M4 Microdebrider	1898200

Contraindications

None known.

Warnings

- W1 Use adequate irrigation from a separate user-provided irrigating source. The use of a bur without irrigation may cause an inordinate amount of heat buildup resulting in thermal injury to tissue.
- W2 Do not use an overheated device, as it may cause thermal injury to the patient or operator.
- W3 When not operating handpiece, eliminate accidental foot control activation. Control energy to and through the handpiece to prevent unintended tissue, bone, or nerve resection.
- W4 Disposable devices are for single-use only. Do not attempt to sterilize disposable blades, burs or irrigation tubing. These disposables are packed sterile and not intended for repeat use. To prevent contamination use only once.
- W5 After each procedure, properly clean all reusable system components.
- W6 Blade and bur accessories are available for resection of soft tissue and bone for surgical procedures. Use of accessories depends on the intended application and patient needs. Sharp-cutting powered accessories induce bleeding and removal of significant tissue and bone.
- W7 Use care in application of moving cutting end to only appropriate anatomical landmarks and the intended surgical site when using blades and burs.
- W8 Ensure adequate visualization when using blades and burs. Discontinue powered application in the event of lack of visualization of the surgical site.
- W9 Use care in application of the moving cutting end to only appropriate anatomical landmarks and the intended surgical site when using accessories. The use of powered reciprocating/rotating instruments may result in vibration-related injury. Use appropriate precautions.
- W10 Always inspect the components before and after use for any damage. If damage is observed, do not use the damaged part until it is repaired or replaced. Damaged parts may deposit metal shavings on surgical site.
- W11 Do not use any parts other than Medtronic components as damage or substandard performance could result.
- W12 Remove and discard accessories following local regulations for proper disposal of contaminated materials.

- W13 Lock disposables properly into handpiece before use.
- W14 Do not modify accessories used with the handpiece. Performance could be diminished with modified accessories.
- W15 Insertion of metal objects in blade or bur tip may cause the blade or bur to break leaving fragments in the wound. The fragments may be difficult to remove causing irritation, inflammation and foreign body response at the surgical site.
- W16 Always keep the cutting tip of the blade or bur away from fingers and loose clothing. Prevent lacerations of user and cross-contamination through compromised glove.
- W17 Sterilize and dry handpiece before storing the system. Decrease likelihood of cross-contamination with timely sterilization.
- W18 Bending or prying may break the blade or bur, causing harm to patient or staff.
- W19 Do not change blades or burs with handpiece running to prevent laceration of user and cross-contamination through compromised glove.
- W20 Discontinue use of curved bur if tip begins to wobble and replace the bur to prevent unintended tissue removal from patient.
- W21 Properly dispose of single-use devices removed from sterile packages. Devices lose sterility upon removal from packaging.
- W22 Do not use accessory if package is opened or damaged. Broken seal offers no protection against cross-contamination.
- W23 All service must be performed by Medtronic qualified personnel only.
- W24 Do not attempt to resharpen used tools. Worn tools should be replaced with new ones frequently to ensure effective cutting and control.
- W25 The Indigo High-speed Otologic drill is designed to work only with Medtronic Xomed burs and Indigo attachments. The use of other burs and attachments may result in sub-standard performance and will void the manufacturer's warranty. Do not use with Midas Rex attachments or burs.
- W26 Do not use chlorine based or corrosive cleaning agents such as bleach, lye, acetone, sodium hypochlorite/bleach, sodium hydroxide, formic acid, or solutions containing glutaraldehyde.
- W27 For handpieces and attachments, the use of a washer-disinfector for cleaning may cause a pre-mature degradation in performance.
- W28 Do not soak/submerge devices.
- W29 Do not use Ultrasonic to clean devices.
- W30 Handpieces and attachments may fail due to extended use and allow a component to detach and fall from the handpieces or attachment, causing patient injury.
- W31 Do not operate the XPS Nexus System in the presence of flammable anesthetics. Avoid potential ignition or explosion of gases.
- W32 Disconnect power to the XPS Nexus Console before cleaning the unit to avoid electrical macro shock.
- W33 Achieve electrical grounding reliability with proper connections. Connect the XPS Nexus console to hospital grade receptacles only.
- W34 When operating the Indigo drill engage the tissue /bone only upon reaching desired speed and stability.
- W35 It is important that the XPS Nexus operator be familiar with this manual, its precautions, procedures and safety issues.
- W36 This medical device complies with EN60601-1-2 safety standard for electromagnetic compatibility, requirements and test. However, if this equipment is operated in the presence of high levels of electromagnetic interference (EMI) or highly sensitive equipment, interference may be encountered and the user should take whatever steps are necessary to eliminate or reduce the source of the interference. Diminished performance may lengthen operating time for anesthetized patient.
- W37 Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this Guide.
- W38 Portable and mobile RF communications equipment can affect Medical Electrical Equipment.
- W39 Use of accessories and cables other than those specified and sold by Medtronic, Inc. may result in increased emissions and decreased immunity of this unit.
- W40 The XPS Nexus should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the XPS Nexus should be observed to verify normal operation in the configuration in which it will be used.
- W41 Repair and/or modification to the XPS Nexus system by anyone other than qualified service personnel may significantly compromise the unit's ability to perform effectively and/or void the equipment warranty.

Precautions

- P1 Blades should be operated in the oscillate mode only. Operating in the forward mode may cause damage to the blade. Burs should be operated in the forward mode only.
- P2 The handpiece should be evaluated prior to each use for suitable operating condition.
- P3 System components should be operated and inspected for damage prior to use. Do not use system if damage is apparent.
- P4 It is recommended that a secondary handpiece be available to minimize any downtime and inconvenience to the surgical staff.
- P5 The StraightShot Magnum II and StraightShot M4 Microdebridors are intended to operate at speeds greater than 6,000 rpm ONLY when used with the XPS High Speed bur line.
- P6 Always ensure that the bur is securely engaged into the handpiece prior to operating the system.

- P7 When operating or testing the Indigo drill handpieces, ensure the bur is properly inserted and locked into the handpiece. Running the drill with the collet unlocked can damage the locking mechanism.
- P8 Do not allow the console to get wet. If liquid enters the console, damage could occur.
- P9 Do not sterilize the console.
- P10 To prevent damage to curved blades and burs, disconnect suction tube prior to changing blade or bur during procedure.
- P11 Do not clean handpieces in ultrasonic cleaner or allow them to be fully immersed in any soaking solution. These procedures may result in damaging the handpieces beyond repair.
- P12 Do not clean the handpiece in an ultrasonic cleaner or cold soak sterilize the handpiece in glutaraldehyde. Do not immerse handpieces, motors, or handpiece cables in any solution, except as detailed in the Cleaning section of this manual. These procedures void the warranty and may damage the handpiece beyond repair.
- P13 Always handle handpieces with care to avoid damage.
- P14 Do not use organic solvents such as acetone or isopropyl alcohol to clean the bur chuck. Use only an enzymatic detergent and distilled water for cleaning after every case or as required.
- P15 Remove the bur from the handpiece before sterilization.
- P16 Sterilize immediately after cleaning. Do not store unless a drying cycle has been performed.
- P17 Always wrap the electrical cord carefully inside the sterilization tray during sterilization to prevent damage to the cord when closing the tray lid.
- P18 Remove the handpiece from the sterilizer immediately after the sterilization cycle is complete.
- P19 Temperatures higher than those stated may be used for handpiece sterilization when necessary to satisfy governmental or health care facility requirements so long as the temperature does not exceed 149°C (300°F).
- P20 Heating above 149°C (300°F) may damage components and will void the warranty.
- P21 Regardless of which type of steam sterilization is used, it is extremely important that the handpiece is rapidly and completely dried before it is stored. Do not operate or store the handpiece unless a drying cycle has been performed. If a vacuum drying cycle is not used following steam sterilization, moisture may be trapped within the handpiece causing corrosion and residue deposits, resulting in premature wear and a reduction in the functional life expectancy of the handpiece. In addition, damage may result if the handpiece is operated with moisture in the electronic connections.
- P22 Store system in a clean dry place.

When the system arrives

1. Verify the contents of the box match the packing slip. If incomplete or damaged, notify Medtronic Customer Service.
2. If container is damaged, or cushioning material shows stress, notify carrier and Medtronic Customer Service. Keep shipping materials for carrier inspection.
3. Save the cartons and packing material. If the instrument is to be shipped the shipping package will provide proper protection.

Attachment assembly

Refer to the "Indigo high-speed otologic drill" and "StraightShot M4 and StraightShot Magnum II" topics for information on how to attach the handpiece.

Irrigation assembly

Refer to the "Indigo high-speed otologic drill" and "StraightShot M4 and StraightShot Magnum II" topics for information on how to connect the irrigation assembly.

Tool operation

Refer to the "Indigo high-speed otologic drill" and "StraightShot M4 and StraightShot Magnum II" topics for information on how to operate the Indigo and StraightShot M4/Magnum II.

System lights and figures

System lights

The XPS Nexus has a green power indicator light that indicates when the system is on.

The XPS Nexus has an amber light located to the right of the green power indicator light. The amber light indicates either the handpiece hall sensor failed, or the connected handpiece has stalled. Refer to the "Troubleshooting" topic for additional information.

System figures

Figure 1. Nexus Console Front

- 1

2

3

4

5
- 6

7

8

9

10

11

12

13
- 14

15

16
- 10

11

12

13

14

15

16
- 1

2

3

4

5

6

7

8

9
- 10

11

12

13

Figure 2. Console Back

- 1

2

3

4
- 1

2

3

4

Figure 3. Nexus Console Connector Panel

- 1

2

3
- 1

2

3

Figure 4. Nexus System Configuration One

- 1

2

3

4

5
- 1


2

3

4

5

Figure 5. Nexus System Configuration Two

- 1 Floor stand Irrigation pole
- 2 Irrigation bag hanger
- 3 Pole clamp

Figure 6. Foot control



Figure 7. Operating Room Setup

- | | |
|---|-------------------------|
| 1 Anesthesia equipment | 6 Microscope |
| 2 XPS Nexus | 7 Surgeon |
| 3 Nursing supplies/Surgical instruments | 8 Electro-Surgical unit |
| 4 Scrub nurse | 9 Anesthesiologist |
| 5 NIM Monitor | 10 Patient |

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Setup the XPS Nexus

1. Verify Operation Room set up (Figure 7 in the "System figures" section). The surgeon may have preferences to the location and visibility.
2. Inspect all components for damage and determine if the system is ready for use.
3. Mount the XPS Nexus and irrigation/coolant bag on the IV pole (Figure 4 or 5 in the "System figures" section).
4. Plug the XPS Nexus into the power source. Position the XPS Nexus so that it does not obstruct the power source for the purpose of disconnecting the main voltage by the power cord.
5. Locate the foot control connection receptical on the connector panel (Figure 3 in the "System figures" section), align the red mark on the connector to the red mark on the console, then insert the connector.
6. Locate the handpiece connection receptical on the connector panel (Figure 3 in the "System figures" section), align the red mark on the connector to the red mark on the console, then insert the connector.
7. Connect the irrigation tubing to the XPS Nexus by doing one of the following:

For the StraightShot Magnum II /StraightShot M4 microdebridors:

- a. Close the shutoff valve on the irrigation tubing.
- b. Locate the pump and lift up the door (Figure 1).
- c. Adjust the tube roller to the up position.
- d. Insert the irrigation tubing with the tubing coming from the IV bag on the left and going to the handpiece on the right.
- e. Snap the pump lock shut.

Note: Ensure the shutoff valve for the tubing is facing to the rear of the console.

- f. Attach the irrigation tubing to the blade.
- g. Open the shutoff valve.

For the Indigo high-speed otologic drill:

- a. Locate the pump and lift up the door (Figure 2).
- b. Insert the pump cartridge.
- c. Snap the pump door shut.

Warning: Ensure the pump cartridge does not crimp the tubing.

8. Turn on the XPS Nexus.
9. Prime the irrigation by holding the handpiece away from the patient and surgical drapes, then press foot control continuously until irrigation flow appears at the tip of the handpiece.

Figure 1. Install StraightShot Irrigation Tubing.



Figure 2. Install Indigo Irrigation Tubing.

Change the default settings on the XPS Nexus

For the StraightShot Magnum II/M4, you can adjust the RPM, switch between oscillate and forward, and adjust the flow rate. The following are the default settings when the StraightShot Magnum II/M4 is attached:

- Oscillation speed = 5000 RPM and the irrigation = 5cc/min.
- If the user changes the mode to Forward, the speed = 12000 RPM and irrigation = 30cc/min.

For the Indigo, the user can adjust the direction (FWD, REV), the RPM, and the flow rate. The following are the default settings when the Indigo is attached:

Forward speed = 50000 RPM and the irrigation = 20 cc/min.

If the user changes the mode to Reverse, the speed = 50000 RPM and irrigation = 20 cc/min.

The following figures describe how to adjust the XPS Nexus.

Figure 1. Adjust the cc/min setting using the Pump panel


Figure 2. Adjust the RPM.


Figure 3. Select oscillation or forward rotation.


Figure 4. Adjust the RPM.



Figure 5. Change the direction of rotation.

Console reuse and cleaning

The XPS Nexus can be cleaned with warm water and mild detergent. Avoid the use of aggressive solvents, bleaches, or sterilization procedures to prevent damage.

Sterility

Do not attempt to resterilize.

Service and maintenance

Contact Medtronic for any service and maintenance inquiries. Do not attempt to fix the XPS Nexus on your own as this is not a user-serviceable device.

There are no maintenance requirements for the XPS Nexus console or foot pedal.


Indigo high-speed otologic drill

Device description

The Indigo drill is a small, compact high-speed, high-torque, reversible electric drill that can be used to dissect bone and biomaterial at variable speeds from 8000 to 60000 rpm. The cable cannot be removed from the drill.

The following instructions for the Indigo drill are in addition to "Set up the XPS Nexus" general assembly instructions. Complete XPS Nexus setup, then continue to the instructions below.

Figure 1. Indigo Alignment

Indigo drill straight attachment assembly

Important: The Indigo High-speed Otologic drill is designed to work only with Medtronic Xomed burs and Indigo attachments. The use of other burs and attachments may result in sub-standard performance and will void the manufacturer's warranty. Do not use with Midas Rex attachments or burs.


1. Verify the alignment marks on the motor collet are in alignment (Figure 1). **Note:** If the marks are misaligned, turn the collet until the marks are aligned.
2. Slide the attachment over the motor collet (Figure 2) so that the alignment mark on the attachment aligns with the alignment mark (unlocked symbol) on the motor collet.
3. Insert the burr with a slight twisting motion until you feel it seat into position.
4. Turn the attachment so the alignment mark aligns with the locked symbol on the motor collet. You will hear two clicks while rotating the attachment.

5. To ensure a secure fit, gently pull the tool.

Important: The Indigo motor will not run correctly unless the attachment is in the locked position.

Warning: Smoke may be generated if attachment is not in the locked position.

Figure 2. Indigo Straight Attachment Assembly

1 Bur (applied part)

Indigo drill angled attachment assembly


Important: The Indigo High-speed Otologic drill is designed to work only with Medtronic Xomed burs and Indigo attachments. The use of other burs and attachments may result in sub-standard performance and will void the manufacturer's warranty. Do not use with Midas Rex attachments or burs.

1. Verify the alignment marks on the motor collet are in alignment (Figure 1).
Note: If the marks are misaligned, turn the collet until the marks are aligned.
2. Slide the attachment over the motor collet (Figure 3) so that the alignment mark on the attachment aligns with the alignment mark (unlocked symbol) on the motor collet.
3. Turn the attachment so the alignment mark aligns with the locked symbol on the motor collet. You will hear two clicks while rotating the attachment.
4. Verify the alignment marks on the tool lock ring align with the unlocked symbol on the attachment.
5. Insert the bur with a slight twisting motion until you feel it seat into position.
6. Turn the lock ring so the alignment mark aligns with the locked symbol on the attachment. You will hear two clicks while rotating the attachment.
7. To ensure a secure fit, gently pull the tool.

Important: The Indigo motor will not run correctly unless the attachment and lock ring is in the locked position.

Warning: Smoke may be generated if attachment is not in the locked position.

Figure 3. Indigo Angled Attachment Assembly


1 Tool Lock Ring

2 Bur (applied part)

Indigo drill irrigation and tubing clip assembly

- 1. Snap or slide the irrigation clip on to the attachment (Figure 4).
Note: The irrigation clip is designed with two nodules that fit in the grooves on the attachment.
- 2. Bend the irrigation tube to a desirable angle and adjust the clip closer to or farther from the bur, as necessary.
- 3. Snap the tubing clip on the handpiece or angled attachment.

Figure 4. Irrigation and Tubing Clip Assembly

Indigo drill cleaning and sterilization instructions

Refer to the "Reprocessing Instructions Indigo High-Speed Otological Drill" and "Reprocessing Instructions Indigo High-Speed Otological Drill Attachment" sections for instructions on cleaning and sterilization.

Indigo drill technical specifications

Indigo Drill 1845000

Size	11.9 cm length x 1.53 cm diameter
Weight	102 g
Speed	60000 rpm forward/reverse
Duty Cycle for Applied Part	For continuous use in operating room temperatures up to 33°C, the Indigo motor is rated for 3 minutes at 60000 rpm, followed by 25 minutes of rest.
	For normal operating room temperatures (typically 20°C), the Indigo motor is rated for continuous cutting indefinitely at 60000 rpm.

StraightShot M4 and StraightShot Magnum II Microdebridors

Device description

The following instructions for the StraightShot M4 and the StraightShot Magnum II are in addition to "Set up the XPS Nexus" general assembly instructions. Complete XPS Nexus setup, then continue to the instructions below.

Figure 1. StraightShot M4 Handpiece



Figure 2. StraightShot Magnum II Handpiece

StraightShot M4 blade or bur assembly

1. Insert the tool aligning the tabs with the notches (Figure 3). Orientate the irrigation barb to the left or right side.

Note: The StraightShot M4 uses a four-tab alignment system.

- For rotating straight blades, orient the irrigation barb at the 3 o'clock position for right-handed surgeons and 9 o'clock for left-handed surgeons.
- For rotating curved blades, orient the irrigation barb at 3 o'clock.
- Adjust the finger wheel with small back-and-forth motions until the blade is engaged.

2. Press the locking collar (Figure 4).

3. Release the locking collar.

Note: If the collar does not return to full out position adjust the finger wheel with small back-and-forth motions until collar pops out.

4. Pull on the blade or the bur to ensure engagement and visually check to make sure the distal tip of the inner blade is in contact with the distal tip of the outer cannula (Figure 5).

Figure 3. M4 Blade or Bur Assembly


Figure 4. M4 Blade or Bur Assembly



Figure 5. M4 Blade or Bur Assembly

StraightShot M4 suction and irrigation tube assembly

1. Attach a suction tube to the suction source and an irrigation tube on the irrigation barb (Figure 6).
2. Secure suction and irrigation in the irrigation groove on the handpiece (Figure 6-7).

Figure 6. M4 Suction and Irrigation Assembly




Figure 7. M4 Suction and Irrigation Assembly

StraightShot Magnum II blade and bur assembly

- 1. Press the collet and insert blade in collet (Figure 8).
- 2. Release the collet (Figure 8).
- 3. Pull on the tool to ensure engagement and check distal tip of inner blade is in contact with the distal tip of the outer cannula.

Figure 8. StraightShot II and Bur Assembly

Microdebrider blade control

Important: If airway blade becomes clogged during use, 1-5cc of irrigant could be aspirated by the patient before you detect the clog.

Note: Periodically submerge blade tip in sterile water, with suction on, to keep blades clear during the procedure.

To rotate the outer blade (Figure 12), use the finger wheel (Figure 1).

Cleaning and sterilization instructions

Refer to the "Reprocessing Instructions Microdebriders" section for instructions on cleaning and sterilization.

StraightShot M4 technical specifications

StraightShot M4 1898200T.


Size	14.3 cm length x 1.8 cm width (1898200T)
Weight	228g 1898200T
Speed	50-5000 rpm oscillate 50-12000 rpm forward
Duty Cycle for Applied Part	Under full load is rated for intermittent operation per the following: Maximum On Time: 55 seconds Minimum Off Time: 30 seconds

StraightShot MAGNUM II technical specifications

StraightShot Magnum II 1897200

Size	17 cm length x 1.6 cm width (1897200/1897201)
Weight	240g 1897200, 1897201 254g 1897200T
Speed	50-5000 rpm oscillate 50-12000 rpm forward
Duty Cycle for Applied Part	Under full load is rated for intermittent operation per the following: Maximum On Time: 55 seconds Minimum Off Time: 30 seconds

Figure 12. Blade Dissection

- 1 Suction flow in through inner blade
Irrigation flow between inner and outer blades
- 2 Outer blade
- 3 Inner blade
- 4 Outer sleeve

Troubleshooting

For any troubleshooting items not corrected by the actions below, contact Customer Service.

XPS Nexus and foot control

Issue	Possible Cause	Action
Pump does not run.	Failed internal components. Moisture in cable conflicts with handpiece recognition.	Contact Customer Service Run a dry cycle when sterilizing.
Little or no irrigation flow.	Tubing Set Improperly seated in pump. Tubing is pinched or kinked. Tubing clamps are restricting flow. Irrigation flow rate setting low. Irrigator obstructed.	Reposition tubing in pump, verify pump lid is fully closed with the fluid flow from left to right. Check tubing at side of pump, see Irrigation/Coolant Pumps Check remaining tubing for pinched or kinked areas, if necessary replace tubing. Set tubing clamps in "open" position. Adjust irrigation flow rate. Replace Irrigator.
Console does not power up.	Power cord not properly connected. No power. Power Inlet Fuses blown. Failed internal components.	Connect power cord. Check power available (i.e. power strip is on, circuit breaker is closed etc.) Replace fuses with 5X20 mm T, 4AH, 250V, delayed fuses. Contact Customer Service.
Console does not power down.	Power switch failure.	Unplug power cord, Contact Customer Service.
Amber light is illuminated as soon as XPS Nexus is powered on.	Magnum II/M4 handpiece has problem.	Replace handpiece.
Amber light is illuminated when handpiece is activated.	Handpiece overload due to strenuous use.	Turn console off then on (retry, if necessary).
Foot control unit pedal does not respond	Connector not fully inserted Internal component failure.	Disconnect and reconnect the FCU cable connector. Try different FCU or console to isolate the problem. Contact Customer Service
Handpiece fails to rotate	Failed foot control. Failed handpiece motor or motor driver.	Contact Customer Service. Contact Customer Service.

MEDTRONIC XOMED blades or burs

Issue	Possible Cause	Action
Appears to be damaged or defective. Tool vibrates excessively, abnormal noise movement.	Damaged or defective. Tool is not firmly seated.	Remove and replace. Microdebridors, pull back locking collet and re-seat the tool. Visao, unlock collar, check/re-seat notch, lock collar.
No suction.	Blade opening is obstructed. Tubing obstructed.	Use stylet to clear blade. Remove blade from surgical site and submerge the blade tip in sterile water with suction connected to the handpiece to evacuate the obstruction. Remove and inspect suction tubing, and if obstructed, remove obstruction, reconnect tubing.
Tool is leaking irrigant.	Tool not seated correctly in collet. Low or no suction.	Check for proper tool insertion by pulling back locking collet, and re-seating tool. See Possible Cause, No Suction.
Tool wobbles in handpiece.	Tool wobbles in handpiece.	Reduce handpiece operating speeds. Use tools that are rated for the console speed selected. Operate handpiece at 50% of full speed for medium, long and X-long burs. Select a new tool.

System requirements and specifications

Console specifications

Functional Standards for Electrical Systems

ANSI/AAMI 60601-1	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance	2005
IEC 60601-1	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance - with amendment 1(2012-07), Corrigendum 1 (2012)	2012
IEC 60601-1-2	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests - Edition 3.0	2007
IEC 60601-1-6	Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral Standard: Usability - Edition 3.0	2010

Physical Dimensions

Console

Size 265 mm W x 185 mm H x 280 mm D

Weight 4.6 kg

Foot control

Size 200 mm W x 62 mm H x 220 mm D

Weight 1.5 kg

Operational Environment

Temperature +10°C to +33°C

Humidity 30% to 75% RH

Barometric Pressure 700 - 1060 hPa

Transport and Storage Environment

Temperature -40°C to +70°C

Humidity Up to 95% RH at 35°C for maximum 3 days

Barometric Pressure 500 to 1060 hPa

Electrical

Input Voltage 100-240 V~ ± 10%

Frequency 50/60 Hz

Power Consumption 180 VA

Internal Fuse 5 x 20 mm, T 4AH, 250 V

Duty Cycle for Applied Maximum On Time 55 Seconds

Part (Magnum II/IV)

Minimum Off Time 30 Seconds for a 30 min of duty cycle. Minimum 1 hour rest before restarting the duty cycle.

Duty Cycle for Applied Maximum On Time 3 minutes at 60000 rpm

Part (Indigo)

Minimum Off Time 25 minutes

Type of protection Class 1
against electrical shock

Degree of protection Type BF
against electrical shock

Power Cord Product Numbers

North America: USA, Barbados, Belize, Bolivia, Canada, Columbia, Ecuador, Venezuela Standard P/N EA600 or 1895820 6 meter P/N EA650 or 189721	United Kingdom, Ireland, Hong Kong, Malaysia, Singapore P/N EA606 or 1895821	Continental Europe: Austria, Belgium, Finland, France, Germany, Greece, Korea, Luxembourg, Netherlands, Norway, Portugal, Spain, Sweden P/N EA602 or 1895822
Argentina P/N EA608	India, South Africa P/N EA607	Switzerland P/N EA601
Australia, New Zealand P/N EA605	Israel P/N EA609	Denmark P/N EA610
	Japan P/N EA603 or 1895823	Italy, Chile P/N EA611

Guidance and manufacturer's declaration – electromagnetic immunity

Part I

Guidance and manufacturer's declaration – electromagnetic immunity – Part I

The XPS Nexus is intended for use in the electromagnetic environment specified below. The customer or the user of the XPS Nexus should assure that it is used in such an environment.

Immunity test	IEC/EN60601-1-2 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the XPS Nexus requires continuous operation during power mains interruptions, it is recommended that the XPS Nexus be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

(Note: UT is the a.c. mains voltage prior to application of the test level.)

Guidance and manufacturer's declaration – electromagnetic emissions

The XPS Nexus is intended for use in the electromagnetic environment specified below. The customer or the user of the XPS Nexus should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The XPS Nexus uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	The XPS Nexus is suitable for use in all establishments, other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings for domestic purpose.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations IEC 61000-3-3	Complies	

Recommended separation distances between portable and mobile RF communications equipment and the XPS Nexus

The XPS Nexus is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the XPS Nexus can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the XPS Nexus as recommended below, according to the maximum output power of the communications equipment.

Rated maximum power of transmitter W	Separation distance according to frequency of transmitter meters		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.10	0.38	0.38	0.73
1.00	1.20	1.20	2.30
10.00	3.80	3.80	7.30
100.00	12.00	12.00	23.0

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Part II

The XPS Nexus is intended for use in the electromagnetic environment specified below.

The customer or the user of the XPS Nexus should assure that it is used in such an environment.

Immunity test	IEC/EN60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the XPS Nexus, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	$d = 1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the XPS Nexus is used exceeds the applicable RF compliance level above, the XPS Nexus should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the XPS Nexus.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Returns

Contact Medtronic Xomed Customer Service at +1 800 874 5797, or your local representative/distributor to obtain a Return Merchandise Authorization number (RMA#) prior to shipping the Product to Medtronic Xomed. Please have the original invoice number or purchase order number available to assist in verifying warranty information. The RMA# should be prominently displayed on the box and included on all paperwork enclosed with the return. All Product returned to Medtronic Xomed should be safely packed in protective wrapping.

Customer must supply the Purchase Order number; the correct shipping and billing address; and a statement of the problem or reason for return.

Customer service information

For further information regarding the use of this product or to report any problems, please contact Medtronic Xomed at +1 800 874 5797, or your local representative/distributor.

Limited warranty

A. This LIMITED WARRANTY provides assurance for the customer who purchases a Medtronic Xomed Product (hereinafter the "Product") that should the Product fail to function to Medtronic Xomed's published specifications during the term of this LIMITED WARRANTY (one year from the date of shipment for new Product, 90 days from date of shipment for refurbished or used Product), Medtronic Xomed will either replace, repair, or issue a credit (adjusted to reflect the age of the Product) for the Product or any portion thereof. This LIMITED WARRANTY is extended only to the buyer purchasing the Product directly from Medtronic Xomed or from its affiliate or its authorized distributor or representative.

B. To qualify for this LIMITED WARRANTY, the following conditions must be met:

(1) The Product must be used on or before its "Use By" or "Use Before" date, if applicable.

(2) The Product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.

(3) Medtronic Xomed must be notified in writing within thirty (30) days following discovery of a defect.

(4) The Product must be returned to Medtronic Xomed within thirty (30) days of Medtronic Xomed receiving notice as provided for in (3) above.

(5) Upon examination of the Product by Medtronic Xomed, Medtronic Xomed shall have determined that: (i) the Product was not repaired or altered by anyone other than Medtronic Xomed or its authorized representative, (ii) the Product was not operated under conditions other than normal use, and (iii) the prescribed periodic maintenance and services have been performed on the Product.

C. This LIMITED WARRANTY is limited to its express terms. THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED WHETHER STATUTORY OR OTHERWISE, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. In no event shall Medtronic Xomed be liable for any consequential, incidental, prospective or other similar damage resulting from a defect, failure, or malfunction of the Product, whether a claim for such damage is based upon the warranty, contract, negligence or otherwise.

D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under legislation governing the sale of consumer goods. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

	Medtronic Reprocessing Instructions Indigo High-Speed Otologic Drill Attachments	Reprocessing Instructions Cleaning and Sterilization instructions are subject to change. Up to date instructions are available online at manuals.medtronic.com.																																															
Warnings and Precautions	- Do not soak/submerge devices. - Do not use ultrasound to clean devices. - Do not use chlorine based or corrosive cleaning agents such as bleach, lye, acetone, sodium hypochlorite/bleach, sodium hydroxide, formic acid, or solutions containing glutaraldehyde. - The use of a washer-disinfector for cleaning may cause a pre-mature degradation in performance. - Allow an adequate cooling period after steam sterilization.																																																
Limitations	Verify functionality prior to re-use.																																																
Point of Use	This product is provided non-sterile and must be cleaned and sterilized before the first use and any reuse.																																																
Containment and Transportation	It is recommended that devices are reprocessed as soon as is practical following use.																																																
Preparation for Decontamination	Follow hospital procedures.																																																
Cleaning: Detergent	A neutral enzymatic to mild alkaline agent (pH 6.0 to 10.5) is preferred. If a washer-disinfector is used, see the instructions supplied with the washer-disinfector machine to select the recommended cleaning agent.																																																
Cleaning: Automated (Do NOT use ultrasonic washer)	- All attachments must be thoroughly rinsed manually with tap water, ensuring all hard to reach areas are rinsed, prior to transfer to automatic washer for processing. - Remove devices from instrument trays before placing into washer baskets. - Orient devices following recommendations of the washer/disinfector manufacturers.																																																
	Recommended Washer Cycle	<table><thead><tr><th>Phase</th><th>Recirculation Time</th><th>Water Temperature</th><th>Detergent Type and Concentration</th></tr></thead><tbody><tr><td>Pre-Wash</td><td>2 minutes</td><td>Cold tap water</td><td>not applicable</td></tr><tr><td>Wash</td><td>5 minutes</td><td>65° (set point)</td><td>Neutral enzymatic to mild alkaline detergent, pH 6.0-10.5</td></tr><tr><td>Rinse</td><td>1 minute</td><td>Hot tap water</td><td>not applicable</td></tr></tbody></table>	Phase	Recirculation Time	Water Temperature	Detergent Type and Concentration	Pre-Wash	2 minutes	Cold tap water	not applicable	Wash	5 minutes	65° (set point)	Neutral enzymatic to mild alkaline detergent, pH 6.0-10.5	Rinse	1 minute	Hot tap water	not applicable																															
Phase	Recirculation Time	Water Temperature	Detergent Type and Concentration																																														
Pre-Wash	2 minutes	Cold tap water	not applicable																																														
Wash	5 minutes	65° (set point)	Neutral enzymatic to mild alkaline detergent, pH 6.0-10.5																																														
Rinse	1 minute	Hot tap water	not applicable																																														
Cleaning: Manual	- Use the detergent, prepared with lukewarm tap water for the cleaning process. - Thoroughly wipe the straight and angled attachments with a cloth dampened with the detergent. - Wipe entire attachment until all gross soil has been removed. - Using a brush wetted with the detergent, brush the attachments to clean the external surfaces and internal connecting surfaces. - For the angled attachment: Wet an appropriately sized cleaning brush (Ø 2.4mm) with the detergent. Insert the brush into the bore at the front of the attachment. Brush the bore to loosen debris trapped inside. Rinse the bore to remove debris. - For the straight attachment: Wet an appropriately sized cleaning brush (Ø 2.4mm) with the detergent. Push the brush through the straight attachment from the rear to front to loosen and remove debris trapped inside. Rinse bore to remove debris. - Place one half of the straight or angled attachment into the detergent. Do not immerse the entire attachment. Gently agitate the attachments in the detergent and actuate any moveable parts. - Place the other half of the attachment into the detergent and repeat. Do not immerse the entire attachment. - Rinse the attachment thoroughly with tap water. Flush both ends to remove detergent. - Thoroughly dry the attachments. An air gun may be used on the straight attachment to blow moisture out from the rear to front. - Using an aerosol spray lubricant (such as Pana Spray), perform the following steps to lubricate attachments: - Holding the can approximately 10-15 cm (3-6 in) away from the attachment, spray all components that move, rotate, or slide with three quick squirts. - Articulate moveable components to ensure proper lubrication. - Remove excess lubricant with a clean cloth.																																																
	Follow hospital procedures.																																																
Packaging	Place devices in instrument tray, and double wrap instrument case with 1-ply polypropylene wrap. In the US, an FDA approved surgical wrap must be used.																																																
Sterilization (Temperatures are minimum required, times are minimum required)	Steam Sterilization <table><thead><tr><th>Cycle</th><th>Gravity</th><th>Pre-Vac</th><th>Pre-Vac (PR/WHO)</th><th>Pre-vac (UK)</th></tr></thead><tbody><tr><td>Temperature</td><td>132°C</td><td>132°C</td><td>134°C</td><td>134°C</td></tr><tr><td>Time</td><td>25 minutes</td><td>4 minutes</td><td>18 minutes</td><td>3 minutes</td></tr><tr><td>Drying</td><td>20 minutes</td><td></td><td></td><td></td></tr></tbody></table> 100% EtO Sterilization Parameters <table><thead><tr><th>Preconditioning</th><th colspan="4">55°C, 70% relative humidity, 30 minutes</th></tr></thead><tbody><tr><td rowspan="5">Sterilization</td><td>Temperature</td><td>55°C</td><td></td><td></td></tr><tr><td>Relative Humidity</td><td>70%</td><td></td><td></td></tr><tr><td>Ethylene oxide concentration</td><td>725 mg/L</td><td></td><td></td></tr><tr><td>Gas exposure time (full-cycle)</td><td>240 minutes</td><td></td><td></td></tr><tr><td>Aeration</td><td>53-57°C, 18 hours</td><td></td><td></td></tr></tbody></table>			Cycle	Gravity	Pre-Vac	Pre-Vac (PR/WHO)	Pre-vac (UK)	Temperature	132°C	132°C	134°C	134°C	Time	25 minutes	4 minutes	18 minutes	3 minutes	Drying	20 minutes				Preconditioning	55°C, 70% relative humidity, 30 minutes				Sterilization	Temperature	55°C			Relative Humidity	70%			Ethylene oxide concentration	725 mg/L			Gas exposure time (full-cycle)	240 minutes			Aeration	53-57°C, 18 hours		
Cycle	Gravity	Pre-Vac	Pre-Vac (PR/WHO)	Pre-vac (UK)																																													
Temperature	132°C	132°C	134°C	134°C																																													
Time	25 minutes	4 minutes	18 minutes	3 minutes																																													
Drying	20 minutes																																																
Preconditioning	55°C, 70% relative humidity, 30 minutes																																																
Sterilization	Temperature	55°C																																															
	Relative Humidity	70%																																															
	Ethylene oxide concentration	725 mg/L																																															
	Gas exposure time (full-cycle)	240 minutes																																															
	Aeration	53-57°C, 18 hours																																															
Maintenance, Inspection and Testing	- Inspect components for any damage before and after each use. If damage is observed do not use the instrument until it is repaired. - After cleaning and sterilization, verify functionality prior to re-use.																																																
Storage	Store with other sterile devices.																																																
Additional Information	None																																																

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing is actually performed using equipment, materials and personnel in the reprocessing facility to achieve the desired result. This normally requires validation and routine monitoring of the process.

Note: All validations performed per current AAMI TIR12, Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers.

Medtronic recommends incineration of devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE/ CJD diagnosis). NHS Estates HTM 2010 Parts 4 & 6: Appendix 2. Items contaminated with TSE Agents and WHO Infection Control Guidelines for Transmissible Spongiform Encephalopathies refers to a TSE decontamination cycle using a steam autoclave at a temperature of 134-137°C for a single cycle of 18 minutes or repeated for a total of six 3-minute cycles.

 Medtronic		Reprocessing Instructions Indigo High-Speed Otologic Drill		Reprocessing Instructions Cleaning and Sterilization Instructions are subject to change. Up to date instructions are available online at manuals.medtronic.com .			
Warnings and Precautions		- Do not soak/submerge devices. - Do not use ultrasound to clean devices. - Do not use chlorine based or corrosive cleaning agents such as bleach, lye, acetone, sodium hypochlorite/bleach, sodium hydroxide, formic acid, or solutions containing glutaraldehyde. - The use of a washer-disinfector for cleaning may cause a pre-mature degradation in performance. - Allow an adequate cooling period after steam sterilization.					
Limitations		Verify functionality prior to re-use.					
Point of Use		This product is provided non-sterile and must be cleaned and sterilized before the first use and any reuse.					
Containment and Transportation		It is recommended that devices are reprocessed as soon as is practical following use.					
Preparation for Decontamination		Follow hospital procedures.					
Cleaning: Detergent		Use a cleaning agent that is suitable for use on aluminum surfaces. A neutral enzymatic to mild alkaline agent (pH 6.0 to 10.5) is preferred. If a washer-disinfector is used, see the instructions supplied with the washer-disinfector machine to select the recommended cleaning agent.					
Cleaning: Automated (Do NOT use ultrasonic washer)		- Remove devices from instrument trays before placing into washer baskets. - Prior to cleaning, cover the drill cable connector end with Handpiece Cleaning Cap, catalog no. 3318520. - Orient devices following recommendations of the washer/disinfector manufacturers. - After completion of the cleaning steps, remove Handpiece Cleaning Cap or other protective components prior to sterilization.					
		Recommended Washer Cycle					
		Phase	Recirculation Time	Water Temperature	Detergent Type and Concentration		
		Pre-Wash	2 minutes	Cold tap water	Not applicable		
		Wash	5 minutes	66° (set point)	Neutral enzymatic to mild alkaline detergent, pH 6.0-10.5		
		Rinse	1 minute	Hot tap water	Not applicable		
Cleaning: Manual		- Prior to cleaning, cover the drill cable connector end with Handpiece Cleaning Cap, catalog no. 3318520. - Wipe all external surface of the motor and cable with a cloth dampened with the detergent prepared with lukewarm tap water. - Brush motor case and collet with a nylon brush dampened with a neutral enzymatic detergent. - Rinse motor thoroughly under running water, collet end pointed down. - Dry collet and motor with lint free towel. - After completion of the cleaning steps, remove Handpiece Cleaning Cap or other protective components prior to sterilization. - Verify that devices are visually clean after manual cleaning.					
Disinfection		Follow hospital procedures.					
Packaging		Place devices in instrument tray, and double wrap instrument case with 1-ply polypropylene wrap. In the US, an FDA approved surgical wrap must be used.					
Sterilization (Temperatures are minimum required, times are minimum required)		Steam Sterilization					
		Cycle	Gravity	Pre-Vac	Pre-Vac (PR/WHO)	Pre-vac (UK)	
		Temperature	132°C	132°C	134°C	134°C	
		Time	25 minutes	4 minutes	18 minutes	3 minutes	
		Drying	20 minutes				
		100% EtO Sterilization Parameters					
		Preconditioning	55°C, 70% relative humidity, 30 minutes				
		Sterilization	Temperature		55°C		
			Relative Humidity		70%		
			Ethylene oxide concentration		725 mg/L		
	Gas exposure time (full-cycle)		240 minutes				
	Aeration		53-57°C, 18 hours				
Maintenance, Inspection and Testing		- Inspect components for any damage before and after each use. If damage is observed do not use the instrument until it is repaired. - After cleaning and sterilization, verify functionality prior to re-use.					
Storage		Store with other sterile devices.					
Additional Information		None					

Notes: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing is actually performed using equipment, materials and personnel in the reprocessing facility to achieve the desired result. This normally requires validation and routine monitoring of the process.

Notes: All validations performed per current AAMI TIR12, Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers.

Medtronic recommends incineration of devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE)/ CJD diagnosis. NHS Estates HTM 2010 Parts 4 & 6 Appendix 2, Items contaminated with TSE Agents and WHO Infection Control Guidelines for Transmissible Spongiform Encephalopathies refers to a TSE decontamination cycle using a steam autoclave at a temperature of 134-137°C for a single cycle of 18 minutes or repeated for a total of six 3-minute cycles.

	Reprocessing Instructions Microdebriders StraightShot M4, StraightShot Magnum II	Reprocessing Instructions Cleaning and Sterilization instructions are subject to change. Up to date instructions are available online at manuals.medtronic.com.		
Warnings and Precautions	• Disconnect the power before cleaning. • Do not fully immerse, or ultrasonically clean, this instrument. • Do not cold soak sterilize this instrument in glutaraldehyde. This will void the warranty. • Do not use organic solvents to clean the bur chuck. • For drill handpiece cleaning, cover handpiece cable connector end with Handpiece Cable Cap, Small, catalog no. 3318510 or Handpiece Cleaning Cap, Universal, catalog no. 3318520. (Note: Use 3318520 for Straightshot M4. Use 3318510 for other handpieces.) • After completion of the cleaning steps, remove Handpiece Cable Cap or other protective components installed prior to cleaning.			
Limitations	After cleaning and sterilization, verify functionality prior to re-use.			
Point of Use	• This product is provided non-sterile and must be cleaned and sterilized before the first use and any reuse. • To remove occasional residual buildup on handpiece cable connector, use a soft brush and isopropyl alcohol.			
Containment and Transportation	It is recommended that instruments are reprocessed as soon as is practical following use.			
Preparation for Decontamination	Remove the bur from the handpiece, otherwise disassembly is not required.			
Cleaning: Automated (Do NOT use ultrasonic washer)	• Remove instruments and equipment from any instrument trays. • Use alkaline (pH ≤ 10.5) or neutral pH detergent. • Do not immerse the handpiece.			
	Auto Wash Pre-Cleaning Instructions - Neutral Detergent			
	• Manually rinse the handpiece, cable, and tubing management clips with running tap water to remove all visible soil from the device. Using a soft bristled brush and a neutral detergent solution, gently and thoroughly clean the entire open lumen of the device. Actuate the collet while rinsing and turn the thumbwheel. Hold the handpiece with collet end down when rinsing and run stream of water through the suction port.			
	• Place devices in the washer with the collet facing downward to facilitate drainage.			
	Auto wash Cycle - Neutral Enzymatic Detergent			
	Phase	Minimum Time	Temperature	Detergent
	Pre-wash 1	02:00	cold tap water	N/A
	Enzyme wash	05:00	< 50°C	neutral pH enzymatic detergent, manufacturer's recommended concentration
	Rinse 1	01:00	hot tap water	N/A
	Thermal rinse	01:00	82°C minimum	N/A
Dry the devices with a clean soft cloth. Visually examine for cleanliness.				
Auto Wash Pre-Cleaning Instructions - Alkaline Detergent				
• Manually rinse the handpiece, cable, and tubing management clips with running tap water to remove all visible soil from the device. Using a soft bristled brush and a neutral detergent solution, gently and thoroughly brush the entire open lumen of the device. Actuate the collet while rinsing and turn the thumbwheel. Hold the handpiece with collet end down when rinsing and run stream of water through the suction port.				
• Place devices in the washer with the collet facing downward to facilitate drainage.				
Auto wash Cycle - Alkaline Detergent				
Phase	Minimum Time	Temperature	Detergent	
Pre-wash 1	02:00	cold tap water	N/A	
Wash 1	05:00	66°C (set point)	alkaline detergent pH ≤ 10.5	
Rinse 1	01:00	hot tap water	N/A	
Thermal rinse	01:00	82°C minimum	N/A	
Dry the devices with a clean soft cloth. Visually examine for cleanliness.				

Medtronic

Reprocessing Instructions

Microdebridors

StraightShot M4, StraightShot Magnum II

 Reprocessing Instructions
 Cleaning and Sterilization Instructions are subject to
 change. Up to date instructions are available online at
manuals.medtronic.com

Cleaning: Manual	• Remove instruments and equipment from instrument tray. • Do not immerse the handpiece. • Wipe the handpiece and cable with neutral detergent applied to a clean, non-abrasive cloth. • Manually rinse the handpiece, cable, and tubing management clips with running tap water to remove all visible soil from the device. Hold the handpiece with collet end down when rinsing and run stream of water through the suction port. • Use an enzymatic detergent solution to loosen and remove collected tissues from the unit. • Using an enzymatic detergent solution and a soft cloth dampened with the detergent, wipe the motor and cable. Gently clean the handpiece, cable, and tubing management clips with a soft bristle brush wetted with the enzymatic detergent. Using a soft bristled brush, gently and thoroughly clean the entire open lumen of the device. Actuate the collet while cleaning and turn the thumbwheel. Use the brush to thoroughly clean the tubing management clips. • Manually rinse the handpiece, cable, and tubing management clips with running tap water to remove all visible soil from the device. • Hold the handpiece with collet end down when rinsing the lumen thoroughly.* * Additional Cleaning Instructions for XPS Straightshot M5/M4/Spine Shaver Microdebrider: • During the normal cleaning cycle, run a gentle stream of warm water into the collet (front end), and into the lock lever of the Straightshot M4/Spine Shaver handpiece. • While warm water is running into the collet, rotate the mechanism for several revolutions (rotate the wheel); and while water is running into the lock lever, actuate the lock lever several times (locking and unlocking). • Shake excess water from the handpiece. • PRECAUTION: Ensure the use of a very gentle stream of warm clean water during this additional cleaning step. • Dry the handpiece and cable with a lint-free towel. Make sure to dry off the electrical connection on the cable ends. • Apply a small amount of silicone spray into the front-end collet and outside of the handpiece. • Sterilize the handpiece immediately after cleaning.																																																																									
Disinfection	Do not cold soak in glutaraldehyde.																																																																									
Packaging	• A standard, sterilization wrap may be used. In the US, an FDA approved surgical wrap must be used. Ensure that the pack is large enough to contain the instrument without stressing the seals. • In sets: Instruments may be loaded into dedicated instrument trays. Wrap trays using appropriate method.																																																																									
Sterilization (Temperatures are minimum required, times are minimum required)	The sterilization parameters given below should be used for devices that are fully disassembled when disassembly is possible. Use basic aseptic technique during post-sterilization assembly to maintain the sterility of the instrument(s). All steam sterilization cycles have been validated in the wrapped configuration and instruments can be sterilized wrapped or unwrapped. <table><tr><th colspan="7">Steam Sterilization</th></tr><tr><th>Cycle</th><th>Gravity</th><th>Gravity*</th><th>Pre-Vac</th><th>Pre-Vac (PR/WHO)</th><th>Pre-vac (UK)</th><th></th></tr><tr><td>Temperature</td><td>132°C</td><td>121°C</td><td>132°C</td><td>134°C</td><td>134°C</td><td></td></tr><tr><td>Time</td><td>15 minutes</td><td>40 minutes</td><td>4 minutes</td><td>18 minutes</td><td>3 minutes</td><td></td></tr><tr><td>Drying</td><td colspan="6">25 minutes, or until visibly dry</td></tr></table> 100% EtO Sterilization Parameters (set points) <table><tr><th>Preconditioning</th><td colspan="6">54°C, 60% relative humidity, 60 minutes</td></tr><tr><th rowspan="5">Sterilization</th><td colspan="3">Temperature</td><td colspan="2">54°C</td><td></td></tr><tr><td colspan="3">Relative Humidity</td><td colspan="2">60%</td><td></td></tr><tr><td colspan="3">Ethylene oxide concentration</td><td colspan="2">725 mg/L</td><td></td></tr><tr><td colspan="3">Gas exposure time (full-cycle)</td><td colspan="2">120 minutes</td><td></td></tr><tr><td colspan="3">Aeration</td><td colspan="2">54°C, 8 hours</td><td></td></tr></table>	Steam Sterilization							Cycle	Gravity	Gravity*	Pre-Vac	Pre-Vac (PR/WHO)	Pre-vac (UK)		Temperature	132°C	121°C	132°C	134°C	134°C		Time	15 minutes	40 minutes	4 minutes	18 minutes	3 minutes		Drying	25 minutes, or until visibly dry						Preconditioning	54°C, 60% relative humidity, 60 minutes						Sterilization	Temperature			54°C			Relative Humidity			60%			Ethylene oxide concentration			725 mg/L			Gas exposure time (full-cycle)			120 minutes			Aeration			54°C, 8 hours		
Steam Sterilization																																																																										
Cycle	Gravity	Gravity*	Pre-Vac	Pre-Vac (PR/WHO)	Pre-vac (UK)																																																																					
Temperature	132°C	121°C	132°C	134°C	134°C																																																																					
Time	15 minutes	40 minutes	4 minutes	18 minutes	3 minutes																																																																					
Drying	25 minutes, or until visibly dry																																																																									
Preconditioning	54°C, 60% relative humidity, 60 minutes																																																																									
Sterilization	Temperature			54°C																																																																						
	Relative Humidity			60%																																																																						
	Ethylene oxide concentration			725 mg/L																																																																						
	Gas exposure time (full-cycle)			120 minutes																																																																						
	Aeration			54°C, 8 hours																																																																						
Maintenance, Inspection and Testing	• Inspect components for any damage before and after each use. If damage is observed do not use the instrument until it is repaired. • After cleaning and sterilization, verify functionality prior to re-use.																																																																									
Storage	It is extremely important that the handpiece be rapidly and completely vacuum dried before storage to prevent corrosion and residue deposits in the bearing and motor.																																																																									
Additional Information	Increase temperatures higher than those stated when necessary to satisfy governmental or health care facility requirements so long as the temperature does not exceed 149° C. Heating above 149° C may damage the handpiece and will void the warranty.																																																																									

* This 40 minutes and 121°C gravity sterilization cycle is not considered by the United States Food and Drug Administration (US FDA) to be a standard sterilization cycle. Users should only use sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization containers) that have been cleared by the US FDA for the selected sterilization cycle specifications (40 minutes and 121°C).

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing is actually performed using equipment, materials and personnel in the reprocessing facility to achieve the desired result. This normally requires validation and routine monitoring of the process.


Note: All validations performed per current AAMI TIR12, Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers.

Medtronic recommends incineration of devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE)/ CJD diagnosis. NHS Estates HTM 2010 Parts 4 & 6; Appendix 2, Items contaminated with TSE Agents and WHO Infection Control Guidelines for Transmissible Spongiform Encephalopathies refers to a TSE decontamination cycle using a steam autoclave at a temperature of 134-137°C for a single cycle of 18 minutes or repeated for a total of six 3-minute cycles.

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A large, stylized signature of the notary, written in dark ink.

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


ATTESTATION

STATE OF FLORIDA
DUVAL COUNTY

On this 19th day of DECEMBER, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following, Annex to IFU(s), issued by Medtronic Xomed, Inc., presented to me by the document's custodian, Alexandra Oliver, and, to the best of my knowledge, that the photocopied document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

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
ATTESTATION

STATE OF FLORIDA
DUVAL COUNTY

On this 19th day of DECEMBER, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following, Annex to IFU(s), issued by Medtronic Xomed, Inc., presented to me by the document's custodian, Alexandra Oliver, and, to the best of my knowledge, that the photocopied document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

Notary:

Judith L. Fisher
Official Notary Signature

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
Goods Receiving
Incoming inspection

Kitting Material

FAIL

QMI Input

OQC2

Pass

Label Preparation (Birth)

MA1

Motor Assembly PSU & Base board Assembly

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apparatus for
disinfection and
sterilization
cleaning

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washer for processing.

Remove devices from instrument trays before placing into washer baskets.

Orient devices following recommendations of the washer/disinfector manufacturers.

Recommended Washer Cycle	Phase	Recirculation Time	Water Temperature	Detergent Type and Concentration
	Pre-Wash	2 minutes	Cold tap water	not applicable
	Wash	5 minutes	66° (set point)	Neutral enzymatic to mild alkaline detergent, pH 6.0-10.5
	Rinse	1 minute	Hot tap water	not applicable

Use the detergent, prepared with lukewarm tap water for the cleaning process.

- Thoroughly wipe the straight and angled attachments with a cloth dampened with the detergent.
- Wipe entire attachment until all gross soil has been removed.
- Using a brush wetted with the detergent, brush the attachments to clean the external surfaces and internal connecting surfaces.
- For the angled attachment: Wet an appropriately sized cleaning brush (Ø 2.4mm) with the detergent. Insert the brush into the bore at the front of the attachment. Brush the bore to loosen debris trapped inside. Rinse the bore to remove debris.
- For the straight attachment: Wet an appropriately sized cleaning brush (Ø 2.4mm) with the detergent. Push the brush through the straight attachment from the rear to front to loosen and remove debris trapped inside. Rinse bore to remove debris.
- Place one half of the straight or angled attachment into the detergent. Do not immerse the entire attachment. Gently agitate the attachments in

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The customer or the user of the XPS Nexus can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the XPS Nexus as recommended below, according to the maximum output power of the communications equipment.

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1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

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Immunity test	IEC/EN60601 1-1-2 test level	Compliance level	Electromagnetic environment - guidance
---------------	------------------------------	------------------	--

Portable and mobile RF communications equipment should be used no closer to any part of the XPS Nexus, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

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Recommended separation distance

Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	$d = 1.2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz

Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:

1 At 80 MHz and 800 MHz, the higher frequency range applies.

2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the XPS Nexus is used exceeds the applicable RF compliance level above, the XPS Nexus should be observed to verify normal operation. If abnormal performance is observed, additional measures may be

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necessary, such as re-orienting or relocating the XPS Nexus
b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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be altered or subjected to misuse, abuse, accident or improper handling.

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D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under legislation governing the sale of consumer goods. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY

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
Handwritten signature of the notary.

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
ATTESTATION


STATE OF FLORIDA
DUVAL COUNTY

On this 19th day of DECEMBER, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following, Technical Information, issued by Medtronic Xomed, Inc., presented to me by the document's custodian, Alexandra Oliver, and, to the best of my knowledge, that the photocopied document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

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



Operating Room Setup


- 1 Anesthesia equipment
 - 2 XPS Nexus
 - 3 Nursing supplies/Surgical instruments
 - 4 Scrub nurse
 - 5 NIM Monitor
 - 6 Microscope
 - 7 Surgeon
 - 8 Electro-Surgical unit
 - 9 Anesthesiologist
 - 10 Patient
- 100-240~, 50/60Hz, 180VA
Nexus Console Front

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



1 Straight shot Magnum II/M4 panel

- 2 Oscillation mode
- 3 RPM Increase/decrease
- 4 RPM selection
- 5 Forward mode
- 6 Indigo panel
- 7 Forward/reverse
- 8 RPM selection
- 9 RPM increase/decrease
- 10 Power light (green)
- 11 Amber light
- 12 Handpiece receptical
- 13 Foot control receptical
- 14 Pump panel
- 15 Flow selection
- 16 Flow increase/decrease

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




Nexus Console Connector Panel


- 1 Handpiece alignment mark (red)
- 2 Foot control receptical

3 Indigo, Straightshot M4 / Magnum II
Чертеж консоли


Стойка «IV Pole» XPS Nexus

3 Indigo, Straightshot M4/Magnum II receptical
Console Drawing

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













Input Current	Input Current: 1,8 Amperes maximum.
Input frequency	50/60 Hz
Input waveform:	Input waveform: Sinusoidal
Display	
The type of sensor	LED lamps under membrane panel
Colors	White, amber and green
Contacting area	Front panel
Physical Dimensions	265 mm W x 185 mm H x 280 mm D Weight 4,6 kg
The SW information	
SW size	60 KB
SW version	V01
SW safety class	B
Release date	April 02, 2015
Duty cycle	Duty Cycle for Applied Part (Magnum II/IV): Maximum On Time 55 Seconds Minimum Off Time 30 Seconds for a 30 min of duty cycle. Minimum 1 hour rest before restarting the duty cycle.

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	Duty Cycle for Applied Part (Indigo): Maximum On Time 3 minutes at 60000 rpm Minimum Off Time 25 minutes
XPS Nexus Foot Control	
The photo	
Press force	0.8 Kilograms
Cord length	3 Meters
Cord diameter	4.5 Millimeters
Physical Dimensions	Size 200 mm W x 62 mm H x 220 mm D Weight 1.5 kg
Power Cord, 6 meter, IEC 320, 115V	
The photo	Refer power cord with the CFN 1895820

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













Medtronic

4.5 - 10.5VDC/7mA


Medtronic Xomed, Inc.,
6743 Southpoint Drive North
Jacksonville, FL 32216 USA

Made in China

REF 1894520

0123

IPX7

M726750C208 A

Medtronic

4.5 - 10.5VDC/7mA

Medtronic Xomed, Inc.,
6743 Southpoint Drive North
Jacksonville, FL 32216 USA

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



	Use with
	Do Not Oil
	Do Not Immerse
	Fuse
IPX1	Protected Against Vertical Water Drop
	Oscillate
IPX7	Protected Against The Effects Of Temporary Immersion In Water
	Precaution: Pinch Hazard. Keep Fingers Clear Of Rollers
ACC	Accessory
	Type BF Applied

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
general and individual residual risks is allowed.

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Step	Action	Acceptance Criteria	Pass/Fail
1.	Connect the console and foot control, power on the console	Power LED indicator light on green	Pass
2.	Connect the Shaver hand piece, check the console can recognize the Shaver hand piece	The Shaver operation zone LEDs shall light on, including rotation direction and speed level	Pass
3.	Check the default setting as below: Mode: Shaver oscillate (OSC) Speed: 5000 rpm Irrigation rate: 5 CC/min	Shaver oscillate LED light on Shaver speed level (5000 rpm) LED light on Irrigation rate (5 cc/min) LED light on	Pass
4.	Set the Shaver speed level at 500 rpm Press the Shaver "+" button, check shaver speed level LED status	Shaver speed increasing from 500, 1500, 3000 to 5000 rpm and stay at maximum level at last. Shaver speed LED light on responding to the shaver speed and stay at maximum level at last	Pass
5.	Press the Shaver "-" button, check shaver speed level LED status	Shaver speed decreasing from 5000, 3000, 1500 to 500 rpm and stay at minimum level at last Shaver speed LED light on responding to the shaver speed and stay at minimum level at last	Pass
6.	Set the irrigation rate at 0 cc/min Press the Irrigation "+" button, check the irrigation rate and LED status	Irrigation rate increasing from 0 (off) start at 5 to 5 cc/min and stay at maximum rate at last Irrigation rate LED light on responding to the selected rate and stay at maximum rate at last	Pass
7.	Press the Irrigation "-" button, check the irrigation rate and LED status	Irrigation rate decreasing from 50 to 5 stop at 0 cc/min and stay at minimum rate at last, when at 0 cc/min, the pump head roller stop	Pass

		Irrigation rate LED light on responding to the selected rate and stay at minimum rate at last	
8.	Press the Shaver forward button, check shaver direction and LED status	Shaver direction is forward Shaver oscillate LED light off Shaver forward LED light on	Pass
9.	Check the default setting as below: Mode: Shaver forward (FWD) Speed: 12000 rpm Irrigation rate: 30 cc/min	Shaver forward LED light on Forward speed level (12000 rpm) LED light on Irrigation rate (30 cc/min) LED light on	Pass
10.	Press the Shaver “-” button, check shaver speed level LED status	Shaver speed decreasing from 12000, 6000, 5000 to 3000 rpm and stay at minimum level at last. Shaver speed LED light on responding to the shaver speed level and stay at minimum level at last.	Pass
11.	Press the Shaver “+” button, check shaver speed level LED status	Shaver speed increasing from 3000, 5000, 6000 to 12000 rpm and stay at maximum level at last. Shaver speed LED light on responding to the shaver speed level and stay at maximum level at last.	Pass
12.	Set the irrigation rate at 0 cc/min Press the Irrigation “+” button, check the irrigation rate and LED status	Irrigation rate increasing from 0 (off) start at 5 to 5 cc/min and stay at maximum rate at last Irrigation rate LED light on responding to the selected rate and stay at maximum rate at last	Pass
13.	Press the Irrigation “-” button, check the irrigation rate and LED status	Irrigation rate decreasing from 50 to 5 stop at 0 cc/min and stay at minimum rate at last, when at 0 cc/min, the pump head roller stop Irrigation rate LED light on responding to the selected rate and stay at minimum rate at last	Pass
14.	Use a clamp to fix the shaver hand piece during it operation and induce overload state, check the console will shut off the shaver hand piece and the	The console shut off the shaver hand piece The amber fault LED indicator lights on	Pass

	amber fault LED will light on		
15.	Turning power off, then on, check the amber fault light condition is reset	The amber LED lights off and the fault light condition is reset	Pass
16.	Use a hall sensor failure simulation fixture to simulate the hall sensor failure state, check the amber fault light condition	The amber LED lights on	Pass
17.	Power off the console, remove the shaver hand piece	N/A	N/A
18.	Power on the console, connect the Drill hand piece, check the console can automatically recognize the Drill hand piece	The Drill operation zone LEDs shall light on, including direction and speed level	Pass
19.	Check the default setting as below: Mode: Drill Forward (FWD) Speed: 50000 rpm Irrigation rate: 20 CC/min	Drill Forward LED light on Drill Speed level (50000 rpm) LED light on Irrigation rate (20 CC/Min) LED light on	Pass
20.	Set drill speed at 20000 rpm Press the Drill "+" button, check drill speed and LED status	Drill speed increasing from 20000, 40000, 50000 to 60000 rpm and stay at maximum level at last Drill Speed LED light on responding to the drill speed level and stay at maximum level at last	Pass
21.	Press the Drill "-" button, check drill speed and LED status	Drill speed decreasing from 60000, 50000, 40000 to 20000 rpm and stay at minimum level at last Drill Speed LED light on responding to the drill speed level and stay at minimum level at last	Pass
22.	Set the irrigation rate at 0 cc/min Press the Irrigation "+" button, check the irrigation rate and LED status	Irrigation rate increasing from 0 (off) start at 5 to 50 cc/min and stay at maximum rate at last Irrigation rate LED light on responding to the selected rate and stay at maximum rate at last	Pass
23.	Press the Irrigation "-" button, check the	Irrigation rate decreasing from 50 to 5 stop at 0 cc/min and stay at minimum rate at last, when at 0	Pass

	irrigation rate and LED status	cc/min, the pump head roller stop Irrigation rate LED light on responding to the selected rate and stay at minimum rate at last	
24.	Shut down the console, then restart, press the Drill Reverse button, check drill direction and LED status	Drill direction is reverse Drill forward LED light off Drill reverse LED light on	Pass
25.	Check the default setting as below: Mode: Drill Reverse (REV) Speed: 50000 rpm Irrigation rate: 20 CC/min	Drill reverse LED light on Drill speed level (50000 rpm) LED light on Irrigation rate (20 CC/Min) LED light on	Pass
26.	Set the drill speed at 60000 rpm Press the Drill “-” button, check drill speed and LED status	Drill speed decreasing from 60000, 50000, 40000 to 20000 rpm and stay at minimum level at last Drill speed LED light on responding to the drill speed level and stay at minimum level at last	Pass
27.	Press the Drill “+” button, check drill speed and LED status	Drill speed increasing from 20000, 40000, 50000 to 60000 rpm and stay at maximum level at last Drill speed LED light on responding to the drill speed level and stay at maximum level at last	Pass
28.	Set the irrigation rate at 0 cc/min Press the Irrigation “+” button, check the irrigation rate and LED status	Irrigation rate increasing from 0 (off) start at 5 to 50 cc/min and stay at maximum rate at last Irrigation rate LED light on responding to the selected rate and stay at maximum rate at last	Pass
29.	Press the Irrigation “-” button, check the irrigation rate and LED status	Irrigation rate decreasing from 50 to 5 stop at 0 cc/min and stay at minimum rate at last, when at 0 cc/min, the pump head roller stop Irrigation rate LED light on responding to the selected rate and stay at minimum rate at last	Pass
30.	Use a clamp to fix the drill hand piece during it operation and induce overload state, check the console will shut off the drill hand piece and the amber fault LED will light on	The console shut off the drill hand piece The amber fault LED indicator lights on	Pass


31.	Turning power off, then on, check the amber fault light condition is reset	The amber LED lights off and the fault light condition is reset	Pass
32.	Complete above steps, make sure during these logic check procedure, all select values shall be indicated by illumination.	All select values were indicated by illumination.	Pass

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Goods Receiving
Incoming inspection
Kitting Material
FAIL
QM Input
OQC2
Pass
Label Preparation (Birth)

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When not operating handpiece, eliminate accidental foot control activation. Control energy to and through the handpiece to prevent unintended tissue, bone, or nerve resection.

Disposable devices are for single-use only. Do not attempt to sterilize disposable blades, burs or irrigation tubing. These disposables are packed sterile and not intended for repeat use. To prevent contamination use only once.

After each procedure, properly clean all reusable system components.

Blade and bur accessories are available for resection of soft tissue and bone for surgical procedures. Use of accessories depends on the intended application and patient needs. Sharp-cutting powered accessories induce bleeding and removal of significant tissue and bone.

Use care in application of moving cutting end to only appropriate anatomical landmarks and the intended surgical site when using blades and burs.

Ensure adequate visualization when using blades and burs. Discontinue powered application in the event of lack of visualization of the surgical site.

Use care in application of the moving cutting end to only appropriate anatomical landmarks and the intended surgical site when using accessories. The use of powered reciprocating/rotating instruments may result in vibration-related

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injury. Use appropriate precautions.

Always inspect the components before and after use for any damage. If damage is observed, do not use the damaged part until it is repaired or replaced. Damaged parts may deposit metal shavings on surgical site.

Do not use any parts other than Medtronic components as damage or substandard performance could result.

Remove and discard accessories following local regulations for proper disposal of contaminated materials.

Lock disposables properly into handpiece before use.

Do not modify accessories used with the handpiece. Performance could be diminished with modified accessories.

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Insertion of metal objects in blade or bur tip may cause the blade or bur to break leaving fragments in the wound. The fragments may be difficult to remove causing irritation, inflammation and foreign body response at the surgical site.

Always keep the cutting tip of the blade or bur away from fingers and loose clothing. Prevent lacerations of user and cross-contamination through compromised glove.

Sterilize and dry handpiece before storing the system. Decrease likelihood of cross-contamination with timely sterilization.

Bending or prying may break the blade or bur, causing harm to patient or staff.

Do not change blades or burs with handpiece running to prevent laceration of user and cross-contamination through compromised glove.

Discontinue use of curved bur if tip begins to wobble and replace the bur to prevent unintended tissue removal from patient.

Properly dispose of single-use devices removed from sterile packages. Devices lose sterility upon removal from

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packaging.

Do not use accessory if package is opened or damaged. Broken seal offers no protection against cross-contamination.

All service must be performed by Medtronic qualified personnel only.

Do not attempt to resharpen used tools. Worn tools should be replaced with new ones frequently to ensure effective cutting and control.

The Indigo High-speed Otolologic drill is designed to work only with Medtronic Xomed burs and Indigo attachments. The use of other burs and attachments may result in sub-standard performance and will void the manufacturer's warranty. Do not use with Midas Rex attachments or burs.

Do not use chlorine based or corrosive cleaning agents such as bleach, lye, acetone, sodium hypochlorite/bleach, sodium hydroxide, formic acid, or solutions containing glutaraldehyde.

For handpieces and attachments, the use of a washer-disinfector for cleaning may cause a pre-mature degradation in performance.

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Do not use Ultrasonic to clean devices.

Handpieces and attachments may fail due to extended use and allow a component to detach and fall from the handpieces or attachment, causing patient injury.

Do not operate the XPS Nexus System in the presence of flammable anesthetics. Avoid potential ignition or explosion of gases. Disconnect power to the XPS Nexus Console before cleaning the unit to avoid electrical macro shock.

Achieve electrical grounding reliability with proper connections. Connect the XPS Nexus console to hospital grade receptacles only.

When operating the Indigo drill engage the tissue /bone only upon reaching desired speed and stability.

It is important that the XPS Nexus operator be familiar with this manual, its precautions, procedures and safety issues.

This medical device complies with EN60601-1-2 safety standard for electromagnetic compatibility, requirements and test. However, if this equipment is operated in the presence of high levels of electromagnetic interference (EMI) or highly sensitive equipment, interference may be encountered and the user should take whatever steps are necessary to eliminate or reduce the source of the interference. Diminished performance

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may lengthen operating time for anesthetized patient.

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this Guide.

Portable and mobile RF communications equipment can affect Medical Electrical Equipment.

Use of accessories and cables other than those specified and sold by Medtronic, Inc. may result in increased emissions and decreased immunity of this unit.

The XPS Nexus should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the XPS Nexus should be observed to verify normal operation in the configuration in which it will be used.

Repair and/or modification to the XPS Nexus system by anyone other than qualified service personnel may significantly compromise the unit's ability to perform effectively and/or void the equipment warranty.

Precautions

Blades should be operated in the oscillate mode only. Operating in the forward mode may cause damage to the blade. Burs should be operated in the forward mode only.

The handpiece should be evaluated prior to each use for suitable operating condition.

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System components should be operated and inspected for damage prior to use. Do not use system if damage is apparent.

It is recommended that a secondary handpiece be available to minimize any downtime and inconvenience to the surgical staff.

The StraightShot Magnum II and StraightShot M4 Microdebridors are intended to operate at speeds greater than 6,000 rpm ONLY when used with the XPS High Speed bur line.

Always ensure that the bur is securely engaged into the handpiece prior to operating the system.

When operating or testing the Indigo drill handpieces, ensure the bur is properly inserted and locked into the handpiece. Running the drill with the collet unlocked can damage the locking mechanism.

Do not allow the console to get wet. If liquid enters the console, damage could occur.

Do not sterilize the console.

To prevent damage to curved blades and burs, disconnect suction tube prior to changing blade or bur during procedure.

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Do not clean handpieces in ultrasonic cleaner or allow them to be fully immersed in any soaking solution. These procedures may result in damaging the handpieces beyond repair.

Do not clean the handpiece in an ultrasonic cleaner or cold soak sterilize the handpiece in glutaraldehyde. Do not immerse handpieces, motors, or handpiece cables in any solution, except as detailed in the Cleaning section of this manual. These procedures void the warranty and may damage the handpiece beyond repair.

Always handle handpieces with care to avoid damage.

Do not use organic solvents such as acetone or isopropyl alcohol to clean the bur chuck. Use only an enzymatic detergent and distilled water for cleaning after every case or as required.

Remove the bur from the handpiece before sterilization.

Sterilize immediately after cleaning. Do not store unless a drying cycle has been performed.

Always wrap the electrical cord carefully inside the sterilization tray during sterilization to prevent damage to the

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2. Inspect all components for damage and determine if the system is ready for use.

3. Mount the XPS Nexus and irrigation/coolant bag on the IV pole.

4. Plug the XPS Nexus into the power source. Position the XPS Nexus so that it does not obstruct the power source for the purpose of disconnecting the main voltage by the power cord.


5. Locate the foot control connection receptical on the connector panel, align the red mark on the connector to the red mark on the console, then insert the connector.

6. Locate the handpiece connection receptical on the connector panel, align the red mark on the connector to the red mark on the console, then insert the connector.

7. Connect the irrigation tubing to the XPS Nexus by doing one of the following:

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
Install StraightShot Irrigation Tubing.

a. Close the shutoff valve on the irrigation tubing.

b. Locate the pump and lift up the door.

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
- c. Adjust the tube roller to the up position.
 - d. Insert the irrigation tubing with the tubing coming from the IV bag on the left and going to the handpiece on the right.
 - e. Snap the pump lock shut.
 - f. Attach the irrigation tubing to the blade.
 - g. Open the shutoff valve.
- For the Indigo high-speed otologic drill:


- a. Locate the pump and lift up the door.
 - b. Insert the pump cartridge.
 - c. Snap the pump door shut.
- Warning: Ensure the pump cartridge does not crimp the tubing.
8. Turn on the XPS Nexus.
 9. Prime the irrigation by holding the handpiece away from the patient and surgical drapes, then press foot control continuously until irrigation flow appears at the tip of the handpiece.

Change the default settings on the XPS Nexus
For the StraightShot Magnum II/M4, you can adjust the RPM,

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switch between oscillate and forward, and adjust the flow rate. The following are the default settings when the StraightShot Magnum II/M4 is attached:

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


- Forward speed = 50000 RPM and the irrigation = 20 cc/min.
- If the user changes the mode to Reverse, the speed = 50000 RPM and irrigation = 20 cc/min.

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Adjust the cc/min setting using the Pump panel

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A neutral enzymatic to mild alkaline agent (pH 6.0 to 10.5) is
preferred. If a washer-disinfector is used, see the instructions
supplied with the washer-disinfector machine to select the

предстерилизационной
очистки /Methods and
apparatus for
disinfection and
sterilization cleaning

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recommended cleaning agent.

All attachments must be thoroughly rinsed manually with tap
water, ensuring all hard to reach areas are rinsed, prior to
transfer to automatic washer for processing.

Remove devices from instrument trays before placing into
washer baskets.

Orient devices following recommendations of the
washer/disinfector manufacturers.

Recommen ded Washer Cycle	Phase	Recirculat ion Time	Water Temperatur e	Detergent Type and Concentra tion
	Pre- Wash	2 minutes	Cold tap water	not applicable
	Wash	5 minutes	66° (set point)	Neutral enzymatic to mild alkaline detergent, pH 6.0- 10.5
	Rinse	1 minute	Hot tap water	not applicable

Use the detergent, prepared with lukewarm tap water for the
cleaning process.

- Thoroughly wipe the straight and angled attachments with a
cloth dampened with the detergent.
- Wipe entire attachment until all gross soil has been removed.
- Using a brush wetted with the detergent, brush the attachments
to clean the external surfaces and internal connecting surfaces.
- For the angled attachment: Wet an appropriately sized cleaning
brush (Ø 2.4mm) with the detergent. Insert the brush into the
bore at the front of the attachment. Brush the bore to loosen
debris trapped inside. Rinse the bore to remove debris.
- For the straight attachment: Wet an appropriately sized

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cleaning brush (Ø 2.4mm) with the detergent. Push the brush through the straight attachment from the rear to front to loosen and remove debris trapped inside. Rinse bore to remove debris.

- Place one half of the straight or angled attachment into the detergent. Do not immerse the entire attachment. Gently agitate the attachments in the detergent and actuate any moveable parts.
- Place the other half of the attachment into the detergent and repeat. Do not immerse the entire attachment.
- Rinse the attachment thoroughly with tap water. Flush both ends to remove detergent.
- Thoroughly dry the attachments. An air gun may be used on the straight attachment to blow moisture out from the rear to front.
- Using an aerosol spray lubricant (such as Pana Spray), perform the following steps to lubricate attachments: Holding the can approximately 10-15 cm (3-6 in) away from the attachment, spray all components that move, rotate, or slide with three quick squirts.

Articulate movable components to ensure proper lubrication. Remove excess lubricant with a clean cloth.

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(4)The Product must be returned to Medtronic Xomed within thirty (30) days of Medtronic Xomed receiving notice as provided for in (3) above.

(5)Upon examination of the Product by Medtronic Xomed, Medtronic Xomed shall have determined that: (i) the Product was not repaired or altered by anyone other than Medtronic Xomed or its authorized representative, (ii) the Product was not operated under conditions other than normal use, and (iii) the prescribed periodic maintenance and services have been performed on the Product.

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D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under legislation governing the sale of consumer goods. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not

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All attachments must be thoroughly rinsed manually with tap water, ensuring all hard to reach areas are rinsed, prior to transfer to automatic washer for processing.

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Recomme nded Washer Cycle	Phase	Recirculat ion Time	Water Temperatur e	Detergent Type and Concentra tion
	Pre- Wash	2 minutes	Cold tap water	not applicable
	Wash	5 minutes	66° (set point)	Neutral enzymatic to mild alkaline detergent, pH 6.0-10.5
	Rinse	1 minute	Hot tap water	not applicable

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- Thoroughly wipe the straight and angled attachments with a cloth dampened with the detergent.
- Wipe entire attachment until all gross soil has been removed.
- Using a brush wetted with the detergent, brush the attachments to clean the external surfaces and internal connecting surfaces.
- For the angled attachment: Wet an appropriately sized cleaning brush (Ø 2.4mm) with the detergent. Insert the brush into the bore at the front of the attachment. Brush the bore to loosen debris trapped inside. Rinse the bore to remove debris.
- For the straight attachment: Wet an appropriately sized cleaning brush (Ø 2.4mm) with the detergent. Push the brush through the straight

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- attachment from the rear to front to loosen and remove debris trapped inside. Rinse bore to remove debris.
- Place one half of the straight or angled attachment into the detergent. Do not immerse the entire attachment. Gently agitate the attachments in the detergent and actuate any moveable parts.
 - Place the other half of the attachment into the detergent and repeat. Do not immerse the entire attachment.
 - Rinse the attachment thoroughly with tap water. Flush both ends to remove detergent.
 - Thoroughly dry the attachments. An air gun may be used on the straight attachment to blow moisture out from the rear to front.
 - Using an aerosol spray lubricant (such as Pana Spray), perform the following steps to lubricate attachments: Holding the can approximately 10-15 cm (3-6 in) away from the attachment, spray all components that move, rotate, or slide with three quick squirts. Articulate movable components to ensure proper lubrication. Remove excess lubricant with a clean cloth.

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Guidance and manufacturer's declaration electromagnetic immunity - Part I			
The XPS Nexus is intended for use in the electromagnetic environment specified below. The customer or the user of the XPS Nexus should assure that it is used in such an environment.			
Immunity test	IEC/EN60601-1-2 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD)	±6 kV contact	±6 kV contact	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
IEC 61000-4-2	±8 kV air	±8 kV air	
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment

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Immunity test	IEC/EN60601 1-1-2 test level	Compliance level	Electromagnetic environment - guidance

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is used exceeds the applicable RF compliance level above, the XPS Nexus should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the XPS Nexus
b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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(5) Upon examination of the Product by Medtronic Xomed, Medtronic Xomed shall have determined that: (i) the Product was not repaired or altered by anyone other than Medtronic Xomed or its authorized representative, (ii) the Product was not operated under conditions other

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than normal use, and (iii) the prescribed periodic maintenance and services have been performed on the Product.

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Medtronic

Medtronic Xomed, Inc.
6743 Southpoint Drive North
Jacksonville, Florida 32216 USA

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Sturdygrade is a registered trademark of Sturdygrade LLC. All rights reserved. Sturdygrade is a registered trademark of Sturdygrade LLC. All rights reserved.