NOTARIZATION

The undersigned:

Remco Bosveld, civil-law notary in Amsterdam, the Netherlands;

Hereby certifies:

The signature on the attached document corresponds with the signature in the Dutch passport with number NYBD0JLF6, valid until 6 maart 2028, of **Aniek Therese Paul Hos**, born in Sittard, the Netherlands, on 16 Augustus 1984.

On the strength of information from the Dutch trade register, supplied today, regarding **Philips Consumer Lifestyle B.V.**, a private company with limited liability (*besloten vennootschap met beperkte aansprakelijkheid*), incorporated under the laws of the Netherlands, having its registered office in Eindhoven, the Netherlands, with address High Tech Campus 52, 5656 AG Eindhoven, the Netherlands, registered in the Dutch trade register (*handelsregister*) under number 17066875 (the 'Company'), I confirm that the said Aniek Therese Paul Hos is a board member (*bestuurder*) of the Company and in this capacity jointly authorised with an other board member to represent the Company.

Furthermore the said Aniek Therese Paul Hos is a proxy holder (*procuratiehouder*) of the Company with the following proxies:

- General, category 1, individually authorized up to an amount of EUR 250,000;
- General, category 3, together with a category 1 authorized up to and amount of EUR 12,500,000; and
- Purchase, category 8, together with the Plantmanager authorized up to an amount of EUR 5,000,000.

The certificate mentioned above constitutes no opinion about the contents of the attached document.

The above-mentioned statements are only based on the English part of the attached document.

Signed in Amsterdam, the Netherlands, on June 17, 2022.

R. Bosveld, civil-law notary

APOSTILLE

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
1. Country: THE NETHERLANDS
This public document
2. has been signed by **mr. R. Bosveld**
3. acting in the capacity of notary at Amsterdam
4. bears the seal/stamp of aforesaid notary

Certified

5. in Amsterdam
6. on 21-06-2022
7. by the registrar of the district court of Amsterdam
8. no. 14532
9. Seal/stamp:
10. Signature:

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
Chapters

Welcome to Lumea	6
Intended use	6
Important safety instructions	7
Contraindications	
Warnings	
Cautions	
Notice	
Clinical study	10
Potential side effects of IPL treatment	11
What is IPL?	13
Lumea overview	14
Attachments	
Treatment schedule	17
Treatment modes	18
Start your treatment	19
Prepare before treatment	
Test your skin	
Perform a treatment	
Clean and store	
After treatment	
Troubleshooting	25
Additional information	27

Welcome to Lumea

Achieve gentle, fast and effective hair reduction with Lumea.

We developed this light-based technology for you, in close cooperation with dermatologists. You can conveniently use Lumea in the comfort of your home. Are you ready to get started?

Scan with your phone camera, QR app, or visit www.philips.com/mylumea

Note: In this user manual, Lumea refers to both the Philips IPL Lumea Series 7000 and Series 8000.

Intended use

Philips Lumea is an over-the-counter device intended for the removal of unwanted hair for a single user. Lumea can be used by persons in the age range of 18-65 years with fair up to medium brown skin tones, and naturally dark blonde, brown; dark brown and black hair. Lumea can be used on the body and the female face (below the cheekbone). Lumea is a handheld device intended for home use.

Is Lumea suitable for you?

Philips Lumea is not designed for everyone. Check the skin tone and hair color tables at the start of this user manual if Lumea is suitable for you.

Suitable body hair colors: Treatment with Lumea is not effective if you have light blond, grey, red or white hairs as hairs with these colors do not absorb enough light.

Suitable skin tones: Dark brown to dark skin tones cannot use Lumea as these skin tones absorb light more easily which could result in pain, skin redness or burns.

Important safety instructions

 **ALWAYS** read and follow these instructions carefully to fully benefit from the use of your Lumea and to minimize the risk of injury, skin reactions and side effects. Make sure to save these instructions for future reference.

Contraindications

DO NOT use Lumea if you have had chemotherapy or radiation therapy in the last 3 months.

DO NOT use Lumea if your sensitivity to pain is reduced after using medication or due to a medical condition.

DO NOT use Lumea if you have a dark brown or darker skin tone.

DO NOT use Lumea if you are pregnant, lactating or breastfeeding.

DO NOT use Lumea if you are under the age of 18 years.

Light sensitivity:

DO NOT use Lumea if your skin is extra sensitive to light, for example if you:

- have a sun allergy;
- are taking certain medication or herbal remedies.

ALWAYS check the instruction leaflet inside your medication package to see if what you're taking can cause increased sensitivity to sunlight.

Medical conditions:

DO NOT use Lumea – or consult your doctor first – if you have:

- epilepsy with flashlight sensitivity;
- an active skin cancer or a history of skin cancer;
- a skin pigmentation disorder such as vitiligo;
- a history of poor wound healing;
- an immunocompromised state caused by a disease such as AIDS or by immunosuppressive drugs;
- a history of collagen disorder, including a history of keloid scar formation.

Warnings

Eye safety - To avoid eye damage

DO NOT use Lumea close to or on the eye, eyelid or eyebrow.

Use the device only below the cheekbones.

Electrical safety - To avoid electrical shock

DO NOT try to open, repair, override or modify Lumea or adapter.

DO NOT use any other attachments or adapter than supplied with Lumea.

DO NOT use Lumea in wet surroundings.

DO NOT let water run into the device.

DO NOT expose Lumea to large temperature changes, especially in a high humidity environment, without giving the device time to acclimate. Moisture can damage the device.

ALWAYS check if the voltage indicated on your Lumea adapter corresponds to the local mains voltage before you connect the equipment.

General safety - To avoid serious injury

DO NOT use Lumea directly over or close to any active implants such as a:

- pacemaker;
- neurostimulator;
- insulin pump.

Check the active implants' manual for the safe distance to other electronic devices.

DO NOT use Lumea if you have reduced physical, sensory or mental capabilities as this could result in incorrect usage.

DO NOT use Lumea, including the adapter cable, within 30 cm of WiFi and/or Bluetooth devices (portable RF communication equipment) or smartphones as it may cause Lumea to malfunction.

ENSURE that the adapter can be easily disconnected from the electrical outlet in case of a malfunction. The adapter is considered a disconnection device.

DO NOT use Lumea next to or stacked with other equipment. This could result in improper operation.

ALWAYS check Lumea before you use it. Do not use the device or adapter if it is damaged. Always replace a damaged part with one of the original type.

DO NOT touch any inner part of a damaged device or adapter.

DO NOT wear the adapter cord around your neck to avoid strangulation.

DO NOT leave Lumea unattended when it is turned on to avoid the risk of fire.

KEEP Lumea out of the reach of children. Children should be supervised to ensure that they do not play with the equipment.

ALWAYS check that the adapter is connected to a properly installed and working wall socket to avoid the risk of fire.

ALWAYS unplug Lumea, including adapter, after use and in case of a power failure to prevent damage to the device.

Cautions

Thermal safety - To avoid pain, skin redness or burns

DO NOT continue the treatment if you experience pain or a skin reaction lasting more than 24 hours. Reduce the light intensity setting or stop using Lumea.

DO NOT use any kind of pencil on the skin you want to treat.

DO NOT use Lumea on skin that is darker than the surrounding treatment area.

These highly pigmented skin areas can lead to excessive heating of the skin.

Darker areas may include:

- dark spots such as moles, warts, large freckles, birthmarks, age spots or dark veins;
- marks added to the skin, such as tattoos or (permanent) makeup.

AVOID exposure to sunlight on the treated areas for at least 48 hours after treatment. After treatment, the skin is more sensitive to UV-light. Always use a SPF50+ sunblock if you do go outside during this period.

DO NOT use Lumea within 4 weeks after sunbathing (natural sun or tanning bed) or within 2 weeks after artificial tanning (creams, sprays or pills).

DO NOT use Lumea on lips, nostrils, ears, nipples, labia minora, perineum, around the anus, penis, or scrotum. These areas are highly sensitive and have thinner and darker skin.

DO NOT use Lumea on the male neck and face. The hairs in these areas have different characteristics.

DO NOT use Lumea with the body or precision attachment on the upper lip or other areas of the face.

DO NOT use Lumea on an active skin condition, as this may worsen the condition or cause additional pain, skin redness or burns. Skin conditions may include:

- Active bruises, rashes, sores, open wounds, scars or sunburns. Wait until the injury is completely healed before resuming treatment.
- Active skin irritations, inflammations or infections. Wait until the skin is completely recovered before resuming treatment.
- Damaged skin after using depilatory (hair removal) creams. Wait at least 2 days before resuming treatment.
- Damaged skin after using skin peeling, (micro)dermabrasion or skin rejuvenation methods. Wait at least 1 week before resuming treatment.

DO NOT subject Lumea to heavy shocks to avoid risk of damage to the lamp.

DO NOT use Lumea if there's lotion, gel or visible hair on your skin or on the device. Only treat skin that's clean, dry and hair-free. Make sure Lumea is clean and dry as well.

DO NOT use Lumea if the light exit window or window frame is damaged.

Check the device for damage before use.

DO NOT flash the same spot multiple times during a treatment. Minimize the overlap of flashes on your skin.

DO NOT treat more often than once every 2 weeks.

DO NOT touch the filter, adapter or attachments directly after use. The parts are hot. Let the parts cool down.

ALWAYS check that your attachment is properly connected to Lumea before use.

DO NOT use Lumea on or near any fabric. Make sure to only treat your skin.

DO NOT use Lumea in direct sunlight. The light can influence the SmartSkin sensor.

Eye safety - To avoid eye discomfort

ALWAYS use Lumea in a well-lit room. Your eyes are more sensitive to stray light in darker conditions.

General safety - To avoid injury

DO NOT share Lumea with other persons for hygienic reasons. Lumea is intended for single person use only.

DO NOT use Lumea over vascular skin disorders (e.g. varicose veins).

Notice

Hair removal by Intense Pulsed Light can cause increased hair growth in some individuals. Body areas at risk include women's jawline, neck and areas above the cheekbone and men's upper back and shoulders. Based upon current data, the highest risk groups for this response are females of Mediterranean, Middle Eastern, and South Asian heritage.

Clinical study

Risks of using Lumea

Out of the 94 subjects who were treated with Lumea in the clinical trial XX experienced side effects: itchiness (X), slight redness (X), swelling (X), burns (X), warm sensation (X), inflammation (X).

Benefits

75 subjects received 4 bi-weekly and 8 treatments once every 4 weeks during a clinical trial. After 4 bi-weekly treatments, an average hair reduction of XX % was achieved. After 8 additional treatments performed once every 4 weeks, an average hair reduction of XX % was achieved.

Potential side effects of IPL treatment


The table below lists side effects that may occur following treatment.

 **DO NOT** continue the treatment if you experience pain or a skin reaction lasting more than 24 hours. Reduce the light intensity setting or stop using Lumea.

Common side effects	What is observed?	How long does it last and how can it be helped?
Itching sensation	Perceived mild itching in the treated area.	These reactions can last up to 24 hours. You can cool the area with a cold pack.
Warm/burning sensation	This reaction is not particularly accompanied with visible skin reactions.	
Skin redness (erythema)	Mild redness without clear edge in the treated area.	
Redness around hair follicles	Mild redness around the hair follicles in the treated area.	These reactions can last up to a week. You can cool the area with a cold pack.
Swelling (edema)	Mild swelling around hair follicles appears as bumps in the treated area, sometimes accompanied with redness.	
Skin dryness	Skin appears rough and flaky in the treatment area.	This reaction can last up to 1 month. You can apply a non-scented moisturizer to the treated area.

Uncommon side effects	What is observed?	How long does it last and how can it be helped?
Burns and blisters	Extensive skin redness, swelling and possibly blisters, accompanied with pain.	These reactions may require medical treatment or medication. We advise you to consult a doctor. Wait with the next treatment until the skin has healed completely and make sure you use a lower light intensity.
Superficial crust	Is a sharply defined brownish and crusty area that can occur on darker skin tones.	
Pigmentary changes/skin decoloration	Skin of the treated area appears darker or lighter than the surrounding skin. Pigmentary changes may occur as a secondary effect of a burn.	Do not treat areas with pigmentary changes/skin discolorations until your skin has regained its normal skin tone. Use a lower light intensity when you start treatment again.
Scarring	A scar may occur as a secondary effect of a burn.	Scars can take longer than a month to heal. Do not treat over scar tissue.
Folliculitis	Appears in the treated areas as swelling and redness around hair follicles, accompanied with itching. It may be combined with pustule formation.	Folliculitis may require medical treatment or medication. We advise you to consult a doctor. Wait with the next treatment until the skin has healed completely.
Skin infection	Appears as inflamed skin. This is a secondary effect which may occur after a skin burn or other injuries.	Skin infection requires medical treatment or medication. We advise you to consult a doctor. Do not treat over skin infections.

What is IPL?

1. IPL stands for Intense Pulsed Light.

Gentle pulses of light are absorbed by the roots of your hair.

2. These pulses of light cause the hair follicles to go into a resting phase.


The hair sheds naturally, and regrowth is prevented.

3. Hair doesn't grow continuously: it grows in cycles.

Not all the hairs in one area are growing at the same time. Since IPL is only effective on hair that's actively growing, it's important to repeat treatments as indicated.


Lumea overview

Not shown:
 - Cleaning cloth
 - Pouch (BRI936, BRI939)

SmartSkin sensor

This determines your skin tone and will help you select a comfortable setting by lighting up one of the light intensity lights on top of your device. It also automatically prevents you from treating areas of your skin that are too dark for safe use.

Ready-to-flash light

The white ready-to-flash light will let you know when you have good skin contact and your Lumea is ready to flash.

If the ready-to-flash light is orange, your skin in this body area is too dark to treat. Treat a lighter colored body area.

Attachments

Your Lumea comes with different attachments you can use on different body areas. Select the attachment that is intended for the area you want to treat.

Body attachment


- largest window size
- a curved-in design to effectively cover and treat areas below the neckline

Facial attachment (all models, except BRI930, BRI931)

- smallest window size for precise treatment below the cheekbones
- small, flat design
- extra integrated filter for comfort

Tip: Place your tongue between your upper lip and your teeth or puff out your cheeks to make the treatment easier.

Precision attachment (BRI937, BRI938, BRI939)

- medium window size
- curved-out design for precise and effective use on bikini and underarm areas

Note: The surface of the attachment that is in contact with your skin is slightly larger than the light exit window. The surface of the attachment should overlap the previous spot to make sure you flash the complete skin area.

Note: Treating over the same area multiple times during the same treatment does not increase hair removal efficacy.

Treatment schedule

Follow the indicated treatment schedule.

24 hours

Test phase

Before treating a new body area, or after recent tanning, do a skin test by following the steps in section 'Test your skin'. Wait 24 hours after the skin test. After the skin test, start the treatment phase.

Treatment phase: 4 bi-weekly treatments

During the treatment phase which lasts 4 treatments, treat yourself every 2 weeks.

Treat yourself within 3 days before or after the planned treatment date during the treatment phase.


Touch-up phase: once every 4 weeks

After the treatment phase, use your Lumea once every 4 weeks to maintain results.

Treat yourself within 4 days before or after the planned treatment date during the touch-up phase.


Note: If you miss a treatment, consult chapter 'Troubleshooting' for guidance. If you deviate from the treatment schedule, you may not achieve the intended results.

Treatment modes

Slide & flash

For larger areas (legs), keep the flash button continuously pressed as you slide Lumea over your skin. Lumea keeps releasing flashes until you release the button or if you lose good skin contact.


Stamp & flash

For small body areas (knees, underarms), place Lumea on your skin and press the flash button. Lift Lumea from your skin and place it beside the spot previously treated. Press the flash button again to release another flash.

Start your treatment

Prepare before treatment

1. Remove hair

Wax, shave or epilate before using Lumea.

Do a treatment before new hair is visible on your skin.

If you choose to wax, wait 24 hours before you use Lumea.


2. Clean and dry skin

No water, oil, cream or body milk should be left on your skin.

3. Check your Lumea for damage

⚠ ALWAYS check Lumea before you use it. Do not use the device or adapter if it is damaged. Always replace a damaged part with one of the original type.

4. Choose the right attachment

Click on the attachment for the area you want to treat.

5. Check your Lumea for dirt


⚠ DO NOT use Lumea if there's any lotion or gel on your skin or on the device. Only treat skin that's clean and dry. Make sure your Lumea is clean and dry as well.

Test your skin

Before your first treatment on each body area or after recent tanning, test your skin's reaction and find your suitable setting for treatment.

1. Follow the steps in section 'Prepare before treatment'.


2. Choose the right attachment

Click on the attachment for the area you want to treat.

3. Plug in. Press the power button

The 5 intensity lights above the power button light up and go out one by one.

4. Place Lumea on your skin

It must sit flush against your skin.

5. Get skin tone setting

One of light intensity indicators will turn white to indicate the proposed setting for you.

The setting is selected automatically after 5 seconds.

You can also press the power button to confirm the setting.

Note: If the ready-to-flash light lights up orange your skin may be too dark to treat, please consult chapter 'Troubleshooting'.

6. Begin skin test

Place Lumea back on your skin.

The ready-to-flash light will light up white.

7. Find your comfortable setting

Press the flash button. If the flash is comfortable, you can increase the setting. If you experience pain, lower the setting. Toggle through the settings by pressing the power button. Change the setting and move Lumea to an adjacent spot. Flash again until you find the highest comfortable setting.

8. Turn off, unplug

Hold the power button for 2 seconds. Unplug the device from the power outlet. Let the device cool down.

9. Clean all parts of your Lumea after use

It is important to clean all parts of your Lumea. Residue like hair, sweat, or dirt can build up from the flashing process and may cause skin reactions if left uncleaned.


10. Wait 24 hours and check your skin for any reactions

If there are no skin reactions after 24 hours, start the treatment with the highest comfortable setting. If the highest comfortable setting resulted in a skin reaction, choose a lower setting that did not result in any skin reactions. If you did not test a lower setting, repeat the skin test.

Perform a treatment

Follow the steps in section 'Prepare before treatment'. Test your skin's reaction and find your comfortable setting before treating a new body area or after recent tanning.

- 1. Choose the right attachment**
Click on the attachment for the area you want to treat.

- 2. Be comfortable**
Make yourself comfortable before you start the treatment.

- 3. Plug in. Press the power button**
The 5 intensity lights above the power button light up and go out one by one.

- 4. Select setting**
Press the power button to toggle through the settings and select your highest comfortable setting that did not result in any skin reactions.

- 5. Begin treatment**
Place Lumea back on your skin.
The ready-to-flash light will light up white.

- 6. Flash. Then move**
 **DO NOT** continue the treatment if you experience pain. Reduce the light intensity setting or stop using Lumea. **DO NOT** flash the same spot multiple times during a treatment. Minimize the overlap of flashes on your skin.

Press the flash button. Move to an adjacent spot and flash again.

- 7. Treat another body area**
If you want to treat another body area make sure you select the right attachment and continue with step 1.

8. **Finished. Turn off and unplug**
Hold the power button for 2 seconds. Unplug your Lumea from the power outlet. Let the device cool down.

9. **Clean all parts of your device after use**
It is important to clean all parts of your Lumea. Residue like hair, sweat, or dirt can build up from the flashing process and may cause skin reactions if left uncleaned.

10. **Next treatment**
Follow the treatment schedule and plan in your next treatment.


Clean and store

- ⚠ ALWAYS** unplug Lumea, including adapter, after use and in case of a power failure to prevent damage to the device.

1. **Cool your Lumea**
After use, unplug your Lumea and let it cool down for 10 minutes.


2. **Click off the attachment**
Gently pull it off.

- ⚠ DO NOT** use the device in wet surroundings.
DO NOT let water run into the device.


3. **Lightly moisten Lumea cleaning cloth**
Dampen the soft Lumea cleaning cloth with a few drops of water.


4. Clean the attachments

Clean each part of the attachment you used, including the metal reflectors.

5. Clean the light exit window

Make sure there's no dirt, residue or moisture left on the surface of the light exit window.

6. Click on an attachment. Then store

Make sure the attachment, including the SmartSkin sensor, is protected from damage. Store all parts in the original packaging or in the provided pouch (BRI936, BRI939) in a dust-free place at the correct environmental conditions.

After treatment

After use, if you want to apply lotions, creams, deodorants, moisturizers or cosmetics to the treated areas, always check the instruction leaflet of the product you want to apply to make sure it is suitable for the condition of your skin at that moment.

Troubleshooting

This chapter summarizes the most common problems you could encounter with the device or if you deviate from your treatment schedule. If you are unable to solve the problem with the information below, visit www.philips.com/support for a list of frequently asked questions or contact the Consumer Care Center in your country as indicated in the included international warranty leaflet.

Common problems while using the device:

Problem

Solution

My Lumea is not flashing / the ready-to-flash light turns orange and starts to blink.

1. Your skin may be too dark. Check the skin tone/hair color chart on the inside cover of this booklet.
2. Make sure your attachment is (properly) connected to the front of your Lumea.
3. Make sure you have full skin contact for at least 5 seconds. In some areas (sensitive and bony area) it might be difficult to make good skin contact. No light should be able to enter between Lumea and your skin. Try repositioning your Lumea.
4. Make sure you press the flash button on the inside handle to treat an area of your skin.

My Lumea produces a strange smell while flashing.

1. Make sure you've removed all hairs in the area you want to treat, before using your Lumea. Burned hairs can cause smells.
2. Make sure your skin is clean and dry before starting a treatment.
3. Make sure your Lumea and its attachments are clean and dry before starting a treatment. Residue like hair, sweat, or dirt can build up from the flashing process and may cause skin reactions if left uncleaned.

The light intensity lights on top of my Lumea are blinking and the ready-to-flash light is orange and blinking.

1. Lumea is in error mode. Press the power button for more than 2 seconds, Lumea should go back to standby. Or unplug Lumea from the mains, it should turn off.
2. Wait 5 minutes. Lumea will automatically go back to standby after 5 minutes without help.

What to do if you deviate from your treatment schedule:

Phase	Problem	Solution
Treatment phase	Missed treatment. <i>I did my treatment more than 3 days after the originally planned date.</i>	Try to stay closer to your planned treatment date. You can continue with the original planned schedule. You may add one additional bi-weekly treatment to your treatment phase. If you deviate from the treatment schedule, you may not achieve the intended results.
	Last bi-weekly treatment. <i>I finished my treatment phase, but I feel I did not reach the claimed results.</i>	You may choose to add an additional bi-weekly treatment to your treatment phase, but claimed results are based on the indicated treatment schedule and individual results vary.
Touch-up phase	Out of treatment window. <i>I did not treat within 4 days of my scheduled treatment date.</i>	Please follow the correct schedule. If you deviate from the treatment schedule, you may not achieve the intended results. You can continue with the original planned schedule or you may choose to restart the treatment phase, but claimed results are based on the indicated treatment schedule and individual results vary.

Additional information

Recycling

- This symbol means that this product shall not be disposed of with normal household waste (2012/19/EU).
- Follow your country's rules for the separate collection of electrical and electronic products. Correct disposal helps prevent negative consequences for the environment and human health.

Electromagnetic compatibility

This product has been tested and found to comply with IEC 60601-1-2 for electromagnetic compatibility (EMC) as Class B according to CISPR 11. Tests levels are listed in the accompanying tables. Only the adapter that is supplied with Lumea should be used with your device.

Lumea is intended for use in the electromagnetic environment specified below. The user of the device should assure that it is used in such an environment to prevent improper operation.

Guidance and manufacturer's declaration – electromagnetic emissions

Emissions	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	Lumea does not use RF energy for its intended function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	Lumea is suitable for use in locations in residential environments and in establishments directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity

Lumea is intended for use in the electromagnetic environment listed here.

The user of the device should assure that it is used in such an environment to prevent improper operation.

Immunity test	Immunity test level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	Floors can be a source of static electricity. When the air is very dry (<20 % RH) the user is advised to use Lumea in a room with concrete or wooden floor (instead of using it in a room with a synthetic floor).
Radiated RF EM fields IEC 61000-4-3	10 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz	±8 kV contact ±15 kV air
Proximity fields from RF wireless Communications Equipment IEC 61000-4-3	385 (27), 450 (28), 710 (9), 745 (9), 780 (9), 810 (28), 870 (28), 930 (28), 1720 (28), 1845 (28), 1970 (28), 2450 (28), 5240 (9), 5500 (9), 5785 (9) MHz (V/m)	Lumea is immune against any RF wireless communication device such as smart phone, Wi-Fi or Bluetooth device at a distance of 30 cm or more. Avoid a closer distance as this may result in improper operation.
RATED power frequency magnetic fields IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	
Proximity magnetic fields IEC 61000-4-39	30 kHz, CW, 8 A/m 134.2 kHz, Pulse modulation 2.1 kHz, 65 A/m 13.56 MHz, Pulse modulation 50 kHz, 7.5 A/m	
Electrical fast transient/burst IEC 61000-4-4	±2 kV 100 kHz repetition Frequency	Mains power quality expected of a typical commercial or residential environment.
Surges IEC 61000-4-5	±0.5 kV, ±1 kV Line-to-line	Mains power quality expected of a typical commercial or residential environment.
Conducted disturbances induced by RF fields IEC 61000-4-6	3 V r.m.s., before modulation is applied 0.15 MHz – 80 MHz 6 V r.m.s., before modulation is applied in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase at 0° 0 % UT; 250/300 cycle	Mains power quality expected of a typical commercial or residential environment. At power interruptions Lumea will switch off.

Note: UT is the AC mains voltage prior to application of the test level.

If in the unlikely event Lumea is locked inoperative in a function by disturbances from other devices beyond the levels stated in the table above, Lumea will shut-off completely and must be restarted by the power button. In that case make sure the distance to the other device causing this disturbance is enlarged.

This device complies with all applicable standards and regulations regarding exposure to electromagnetic fields.

Technical specifications

Electrical specification	The combination of adapter and Lumea is considered: Class II ME EQUIPMENT
Attachment	Type BF applied part The applied part (attachment) can become hot during use <60 °C. If this occurs, we recommend letting the device cool down and limit the contact duration <1 minute.
Mode of operation	Continuous operation
Intended operator	The patient is the intended operator
IP rating adapter	IP22
Light output	
Risk group classification	Exempt
Wavelength range of emitted optical radiation	Body and precision attachment: 520 nm – 1200 nm Facial attachment: 590 nm – 1200 nm
Spot size	Body attachment: 4.1 cm ² (32.7 mm x 12.6 mm) Precision attachment: 3 cm ² (30.0 mm x 10.0 mm) Facial attachment: 2 cm ² (25.0 mm x 8.0 mm)
Optical exposure	Setting 1: 2.5 +/- 0.6 J/cm ² Setting 2: 3.1 +/- 0.6 J/cm ² Setting 3: 3.8 +/- 0.6 J/cm ² Setting 4: 4.4 +/- 0.6 J/cm ² Setting 5: 5.0 +/- 0.6 J/cm ²
Maximum optical energy	Body attachment: 23.0 J Precision attachment: 16.8 J Facial attachment: 11.2 J
Optical homogeneity	Max. +/-30 % deviation from average optical exposure in treatment area
Pulse sequence	Single pulse, for all settings
Pulse interval	1.0 - 2.4 s, for all settings
Pulse width	1.1 +/-0.4 ms, (FWHM) at the highest setting
Transport and storage conditions	Lumea remains operational in NORMAL USE within its specification after transport or storage in the following environmental range: - a temperature range of -25 °C to 70 °C; and - a relative humidity up to 95 % RH, noncondensing; and - an atmospheric pressure range of 700 hPa to 1060 hPa. If the device has been stored in a hot or cold environment, place it in an environment with a temperature of 20 °C for 30 minutes to let it reach a temperature within the usage conditions (5 °C to 35 °C) before you use it.
Operating conditions	Lumea complies with its specifications when operated in NORMAL USE under the following environmental operating conditions: - a temperature range of 5 °C to 35 °C; - a relative humidity range of 15 % to 90 % RH, noncondensing; and - an atmospheric pressure range of 700 hPa to 1060 hPa.

Adapter	
Adapter input	100-240 V ~ 50/60 Hz 1.5 A
Adapter output	19.5 V $\approx$ 3.33 A 65.0 W
Average active efficiency	≥ 88.0 %
Efficiency at low load (10 %)	≥ 79.0 %
No-load power consumption	< 0.21 W
Device weight	Adapter ~ 480 g Handpiece ~ 430 g
Device size	Adapter ~ 10.8 x 6.5 x 6.8 cm Handpiece ~ 20.0 x 14.5 x 9.5 cm
Service life	5 years
Materials	
Housing	PC (Polycarbonate), ABS (Acrylonitril-butadiene-styrene), PA (Polyamide)
Cable	TPE (Thermoplastic elastomer)

Explanation of symbols

The warning signs and symbols are essential to ensure that you use this product safely and correctly and to protect you and others from injury. Below you find the meaning of the warning signs and symbols on the label and in the user manual.

IEC 60417-5333 Type BF applied part
Indicates that the part of the appliance that comes into physical contact with the user (also known as the applied part) is of type BF (Body Floating) according to IEC 60601-1. The applied part is the attachment.

ISO 7000-2620 Humidity limitation
Symbol to identify the acceptable upper and lower limits of relative humidity for transport and storage to which the device can be safely exposed.

ISO 7000-2621 Atmospheric pressure limitation
Symbol to identify the acceptable upper and lower limits of atmospheric pressure for transport and storage to which the device can be safely exposed.

ISO 7000-0632 Temperature limit
Symbol to identify the maximum and minimum temperature limits at which the device shall be stored, transported or used.

ISO 7000-0626 Keep dry
Symbol to indicate that the device shall be kept dry.

IP22

Indicates the protection against ingress of solid foreign objects ≥ 12.5 mm and against harmful effects due to the ingress of dripping water when tilted at 15 degrees.

IEC 60417-5172 Class II equipment
Symbol for 'Class II Equipment'. The adapter is double-insulated.

ISO 7000-3082 Manufacturer
This symbol indicates the legal manufacturer.

ISO 7000-2497 Country of manufacture
This symbol indicates the country at which the product was manufactured and the manufacturing date. To identify the country of manufacture the 'CC' shall be replaced by the two letter country code.

ISO 7000-2492 Batch number
This symbol indicates the manufacturer's batch number to identify a device.

ISO 7000-2498 Serial number
This symbol indicates the manufacturer's serial number to identify a device.

IEC 60417-5032 Alternating current
This symbol indicates alternating current.

IEC 60417-5031 Direct current
This symbol indicates direct current.

ISO 7010-M002 Refer to instruction manual/booklet
This symbol indicates that you have to read the user manual carefully before you use the device. Save the user manual for future reference. (The background color is blue, the color of the symbol is white.)

ISO 7010-W017 Warning; Hot surface
This symbol indicates that you should take care to avoid coming into contact with a hot surface. (The background color is yellow, the color of the symbol is black.)

This symbol indicates the maximum temperature at which the transformer may be operated continuously under normal conditions of use.

IEC 60417-5957 For indoor use only
This symbol means that the device is designed for indoor use.

IEC 60417-5947 Safety isolating transformer, short-circuit proof
This symbol means that the transformer is designed to be short-circuit proof.

IEC 60417-6190 Power supply unit, switch mode
This symbol means that the transformer is considered a Switch Mode Power Supply (SMPS) unit.

IEC 60417-6343 Maximum altitude
To indicate that the appliance is intended to be usable up to the maximum altitude 3000 m.

ISO 7000-3706 Single patient multiple use
This symbol indicates that the device may be used multiple times on a single user.

This symbol indicates that an App is available for use with this device. App Store is a service mark of Apple Inc.

This symbol indicates that an App is available for use with this device. Google Play and the Google Play logo are trademarks of Google LLC.

ISO 7010-W001 General warning sign
This symbol means that the user should consult the instructions for use for important safety information such as contraindications and warnings that cannot, for a variety of reasons, be presented on the device itself. (The background color is yellow, the color of the symbol is black.)

ISO 7000-0434A Caution
This symbol means that the user should consult the instructions for use for important safety information such as cautions that cannot, for a variety of reasons, be presented on the device itself.

A recycling symbol to identify the material of which an item is made, to facilitate recycling or other reprocessing. The symbol may include a number and/or abbreviation.

ISO 7000-1135
This symbol means an object is capable of being recycled – not that the object has been recycled or will be accepted in all recycling collection systems.

Indicates the appliance is a medical device.

Indicates the Unique Device Identifier that is used to identify a specific product.

IEC 60417-6352
This symbol indicates that the adapter should not be used if the plug pins are damaged.

Certification mark that shows that this product conforms to all technical regulations of the Eurasian Customs Union.

If you need information or support, please visit www.philips.com/support or read the international warranty leaflet.

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All rights reserved.