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
Instruction for Use Guidewires

DESCRIPTION

The PREFACE braided guiding sheath has a large non-tapered lumen that allows for intravascular passage of catheters and infusion of heparinized normal saline. The sheath, with sideport tubing and valve, has a pre-shaped distal tip section with radio-paque brite tip. The guiding sheath is pre-shaped to facilitate catheter positioning.

Color coding below indicates the French size of the largest intravascular device that will pass through the braided guiding sheath.

French	Color
8	Blue

Dilator

A dilator is provided to facilitate the percutaneous entry of the braided guiding sheath by forming an atraumatic transition from the skin through the subcutaneous tissue to the vessel.

Guidewire

A guidewire is provided for maintaining access to the puncture site upon removal of the sheath.

INDICATIONS

The intended use of the percutaneous **PREFACE** Braided Guiding Sheath is for the introduction of intravascular electrophysiology catheters into any cardiac chamber.

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None known.

COMPLICATIONS

Procedures requiring percutaneous catheter and/or sheath introduction should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during or after the procedure.

Possible complications include, but are not limited to the following:

- air embolism
- hematoma at the puncture site
- infection
- perforation of the vessel wall.

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Instructions for Use

A. For Use as a Primary Sheath

1. Inspect the braided guiding sheath and accessories for defects. Do not use any defective devices.

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B. For Use as a Transseptal Sheath

1. Inspect the braided guiding sheath and accessories for defects. Do not use any defective devices.
2. Before proceeding, verify that the guiding sheath is of proper size to allow passage of desired catheter.
3. Flush the dilator, guiding sheath and infusion sideport with stopcock with heparinized saline.
4. Slowly insert the entire length of dilator through the valve and into the guiding sheath, snapping it into place.
5. Follow standard practice for vessel puncture, guidewire insertion, and guiding sheath use.
6. Using a gentle twisting motion, advance the assembly over the guidewire and through the tissue and into the vessel. Continue to advance the assembly until the distal tip is positioned in the high superior vena cava/innominate vein area.
7. Follow standard practice of using a continuous drip of anticoagulant fluid under pressure through the sideport when the introducer is in the vessel.
8. Remove guidewire and slowly insert transseptal needle through the guiding sheath/dilator, until the needle tip is seen fluoroscopically to be approximately 5 mm within the dilator tip.
9. Follow standard practice for transseptal puncture and guiding sheath use.
10. Once the guiding sheath is securely located in the left atrium, withdraw the needle 1-2 centimeters into the dilator so that the needle tip is not exposed, and

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slowly withdraw the dilator/needle assembly.

11. Follow manufacturer's recommendations for the catheter or device being introduced via the introducer.
12. When removing the braided guiding sheath, re-insert the dilator and wire fully into the guiding sheath to aid in straightening the tip portion. Then, remove the dilator, wire, and guiding sheath as a unit.

Caution: Withdraw the needle into the dilator prior to removal as the needle tip may cause damage to the sheath, arterial or venous structures.

C. For Use in Exchange of a Previous Sheath Refer to **PERRY** Exchange Dilator Instructions for Use.

1. Insure that the initial sheath in place is flushed and clear, observing standard precautions to prevent air embolism.
2. Inspect braided guiding sheath and accessories for defects. Do not use any defective devices.
3. Before proceeding, verify that the guiding sheath is of proper size to allow passage of desired catheter.
4. Flush the dilator(s), guiding sheath and infusion sideport with stopcock with heparinized saline.
5. Insert the long guidewire into the **PERRY** Exchange Dilator only.
6. Introduce the wire and dilator into the initial sheath, keeping the wire well beyond the exchange dilator tip to prevent complications.

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Caution: Upon withdrawal of dilator, catheter, or other inserted device, aspiration through the sideport is recommended

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to remove any fibrin accumulation.

14. Follow manufacturer's recommendations for the catheter or device being introduced via the introducer.

15. When removing the braided guiding sheath, re-insert the standard dilator and wire fully into the guiding sheath to aid in straightening the tip portion. Then remove the dilator, wire, and guiding sheath as a unit.

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WARNINGS

Discard the guiding sheath after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. Sheaths are extremely difficult to clean after exposure to biological materials and may cause adverse patient events if reused. Accordingly, Biosense Webster will not be responsible for any direct or consequential damages or expenses resulting from reuse of the sheath.

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Injection of contrast media is not recommended.

PRECAUTIONS

- Store in a cool, dark, dry place.
- Do not use open or damaged packages.
- Do not leave a sheath in place for extended periods of time without a catheter or a dilator to support the cannula wall.
- Do not alter this device.
- Use prior to the "Use By" date.
- Do not attempt to advance or withdraw the guidewire and/or catheter sheath introducer if resistance is felt. Stop the procedure and use fluoroscopy to determine the cause.

Failure to do so may damage the sheath or catheter, or cause the gasket cap to detach.

- If significant resistance is encountered during introduction, use a 10F or larger introducer.
- Exposure to temperatures above 54°C (130°F) may damage the sheath.
- Do not expose to organic solvents.
- Inspect the braided guiding sheath before use to verify that its size, shape and condition are suitable for the specific procedure.
- Extreme care must be taken to avoid damage to the vasculature through

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which the guiding sheath passes. The guiding sheath may occlude smaller vessels. Care must be taken to avoid complete blood flow blockage.

- Large internal diameter guiding sheaths require less force on the syringe during injection.
- Should the guiding sheath become severely kinked, withdraw the entire system (guiding sheath, guidewire, if present, and catheter sheath introducer).
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

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Disclaimer of Warranty and Limitation of Remedy

There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Biosense Webster, Inc. product(s) described in this publication. Under no circumstances shall Biosense Webster, Inc. be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind Biosense Webster, Inc. to any representation or warranty except as specifically set forth herein.

Descriptions or specifications in Biosense Webster, Inc. printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Biosense Webster, Inc. will not be responsible for any direct, incidental, or consequential damages resulting from re-use of the product.

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Medical device is stored in a cool, dark, dry place at a room temperature, at a place protected from direct sunlight at temperature from +5 to +25 °C and air relative humidity from 20% to 80% (without condensation).

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STERILIZATION

Sterilized with ethylene oxide gas. Nonpyrogenic. Radiopaque. For single use only. Do not resterilize.

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If the products have come into contact with blood and / or body fluids, disposal and / or destruction according to SanPiN 2.1.7.2790-10 (class B).

Unused products (which had no contact with blood and / or body fluids), including expired, disposal and / or destruction, according to SanPiN 2.1.7.2790-10 (Class A).

Recycle components, or dispose of the product and its residual elements or waste items in accordance with local laws and regulations

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Instruction for Use Cables

Federal (USA) law restricts this device to sale by or on the order of a physician.

DESCRIPTION

These cables provide a means to interface a BIOSENSE WEBSTER catheter to the appropriate equipment. The cables may be reused subject to the cleaning and sterilization restrictions herein.

Codes 39E59M, 39F12M and 39E58M can also be used as an extension cable to facilitate catheter connection to equipment that is located out of the immediate area of use.

INDICATIONS FOR USE

This cable provides a means to interface a Biosense Webster electrophysiology catheter to the appropriate equipment. This cable may be reused subject to the cleaning and sterilization restrictions herein.

CONTRAINDICATIONS

There are no known contraindications to this cable.

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INSTRUCTIONS FOR USE

For information on how to connect this cable, refer to the applicable catheter and equipment instructions for use.

The information below is applicable for codes 39E59M, 39F12M and 39E58M and applies only to the BIOSENSE WEBSTER interface cable used with the BIOSENSE WEBSTER Diagnostic/Ablation temperature sensing (thermocouple) catheter and a compatible RF generator (refer to the specifications below).

1. Use this interface cable to connect the BIOSENSE WEBSTER Diagnostic/Ablation Deflectable Tip catheter to a compatible RF generator.
2. Refer to the RF generator manual for detailed generator operating instructions for RF catheter ablation

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WARNINGS

- Do not immerse the cable connectors in fluids.
- Patient or operator injury can result from improper handling of the cable. To prevent damage, store the cable in a clean and secure area.
- If used in the presence of electrical equipment, noise could be induced into the cable.
- The cable should be tested for leakage (lead-to-lead) and resistance (end-to-end) before reuse (leakage should not exceed 1.0 μA at 250 V and resistance should not exceed 10 Ω for all codes except D134401 where it should not exceed 15 ohms.)
- If a break occurs in the cable wire or the cable becomes otherwise electrically discontinuous, arcing may occur in the patient-return or active circuit and may burn the patient or create a fire

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- Do not expose cable to organic solvents.

Do not use radiation sterilization, as this method will damage these devices.

- Use the cable prior to the "Use By" date on the package label.

- If conventional line powered equipment is used near the patient, follow electromechanical device (IEC 601) and hospital standards to prevent extraneous electrical currents from reaching the electrophysiological equipment, catheter and heart. Such currents may cause lethal arrhythmias.

- The cable should be visually inspected for any damage (i.e. kinks, nicks, crushed or elongated sections) prior to use or reuse. Do not attempt to repair any damages. In case of doubt, discard the cable and do not reuse

- Automated cleaning of the cable is not recommended, and may damage the cable.

- Failure to properly clean cables could lead to inadequate resterilization

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PRECAUTIONS

- Ensure each cable end is connected to the proper equipment. The "system" end of the cable states "To System" and must be connected to the System. The "catheter" end of the cable states the catalog number and must be connected to the catheter. For information on how to connect this cable, refer to

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the applicable catheter and equipment instructions for use.

- Use the cable prior to the “Use By” date on the package label.

- The cable should be visually inspected for any damage (i.e. cuts, kinks, nicks, crushed or elongated sections) prior to use or reuse.

- Do not attempt to repair any damages. In case of doubt, discard the cable and do not reuse.

- To prevent damage store product in a clean and secure area.

The cables are designed for a specific amount of sterilization cycles, as mentioned in “Sterilization” section. Exceeding the amount of uses may impair the proper functioning of the cable. The cable should be discarded after the maximum amount of sterilization cycles have been reached.

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The manufacturer warrants conformance of this medical device to the declared characteristics within shelf life (sterility term) provided that all requirements of the manufacturer’s use instructions are met.

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There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Biosense Webster, Inc. product(s) described in

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this publication. Under no circumstances shall Biosense Webster, Inc. or its affiliated companies, be liable for any special, direct, incidental, consequential or other damages other than as expressly provided by specific law.

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Medical device is stored in a cool, dark, dry place at a room temperature, at a place protected from direct sunlight at temperature from +5 to + 25 °C

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and air relative humidity from 20% to 80% (without condensation).

The medical device is transported in cool, dry, dark place, at temperature from - 40 °C to +54 °C and air relative humidity from 20% to 80% (without condensation).

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STERILIZATION

Cables with the following product codes; 39E59M, 39F12M and 39E58M may be provided sterile or non-sterile. If non-sterile, it may be sterilized using a standard ethylene oxide (EtO) gas process. Do not use either radiation or steam sterilization as these will damage the cable.

The cables may be re-sterilized for re-use. The cable and connector surfaces should be wiped clean with soap and water. The connector should not be immersed in fluids. The recommended maximum number of sterilization cycles is ten (10). Lead-to-lead leakage and end-to-end resistance should be checked after each use. Leakage should not exceed 1.0 microamp at 250 volts and resistance should not exceed 10 ohms. Exceeding ten uses may impair the proper functioning of the cable.

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packaged sterile. The cables may be re-processed and sterilized using steam autoclave, STERRAD or ethylene oxide (EtO) gas. After use, clean the cable thoroughly according to hospital guidelines. The cable can withstand cleaning with soap solutions or alcohol and may be left in contact with cold sterilization solutions for up to two hours.

The remaining cables are provided in a sterile state and may be resterilized for up to 20 uses per the instructions below. These cables may be re-sterilized for re-use using the sterilization modality and parameters identified in Table 1. The cable and connector surfaces should be wiped clean per hospital guidance. The cable connector ends should not be immersed in fluids. Lead-to-lead leakage and end-to-end resistance should be checked prior to each use. Leakage should not exceed 1.0 microamp at 250 volts. Resistance should not exceed 10 ohms, except for code D134401 where it cannot exceed 15 ohms. Visually inspect all cable surfaces in adequate lighting to ensure complete removal of soil/residues after cleaning. The recommended maximum number of sterilization cycles is twenty (20)*. Exceeding twenty sterilization cycles may impair the proper functioning of the cable.

*Sterilization data on file

For the codes D128622, D128709, D128624, CY1212CT, CB3434CT, CY1210CT, CR3434CT, CG2025CT, CR3425CT, CB3412CT and CB3410CT at the point of use, remove excess contaminants with a disposable cloth or paper wipe. Cables should be resterilized as soon as is practical after use.

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It is recommended that the cables are packaged in a breathable pouch prior to resterilization. Cables should be coiled into a loop when placing into the pouch. Use a pouch large enough to prevent

stress to the pouch seams or excessive bending of the cable.

The cable may be resterilized using the methods specified in Table 1. Sterilization cycles should be performed using techniques and parameters compliant with the equipment manufacturer's

instructions and facility protocols.

The resterilization parameters provided in Table 2 have been validated by the manufacturer.

- These cables meet the requirements of ANSI/AAMI EC53:19982 and ANSI/AAMI EC53A:1998/(R) 20082.
- Recommended parameters are validated for device functionality. Efficacy of cleaning and sterilization is the responsibility of the user.
- Methods described here are for reference only and no guarantee is made that effective sterilization will result.
- Chamber dimensions, temperature distribution profile and the configuration of products placed in the chamber are variables that must be considered in verifying the effectiveness of a specific sterilization cycle.

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Sterilize by procedures validated for effectiveness, applying suitable biological controls

1. Temperature should not exceed 140°C.
2. At ambient temperatures and % relative humidity conditions.
3. Refer to STERRAD Operator's Manual for operating instructions.

Efficacy of cleaning and sterilization is the responsibility of the user. Methods described here are for reference only and no guarantee is made that effective sterilization will result.

Do not reuse the cable after ten (10) uses. The cable is designed for a maximum of ten sterilization cycles. Exceeding ten uses may impair the proper functioning of the cable. The cable should be discarded after the ten resterilization cycles have been reached. Technomed Europe will not be responsible for any direct or consequential damages or expenses that result from reusing the cable beyond the indicated reuse limitations.

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expired, disposal and / or destruction, ac-
cording to SanPiN 2.1.7.2790-10 (Class
A).

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Recycle components, or dispose of the product
and its residual elements or waste items in ac-
cordance with local laws and regulations.

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These medical products - disposable.
Maintenance and repair are not subject.

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The codes 39F54R; 39F55R; 39F40R;
39F41R; 39F43R; 39F44R produced by
Technomed Europ for Biosense Webster Inc.

Biosense Webster, Inc.

33 Technology Drive Irvine,
California 92618USA

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Instruction for Use External Reference-Patches

- NON-STERILE.
- For single use only.
- Do not use if package is open or damaged.
- Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

DESCRIPTION

The REFSTAR RMT External Reference Patch is an integral part of a nonfluoroscopic catheter tip location and electrogram capture technology called CARTO RMT System. When used with the CARTO RMT System, the intracardiac location of the NaviStar RMT Mapping Catheter tip is compared to the location of this electrode-patch (referred to as the Reference Patch). This reference device consists of a sensor embedded in an adhesive patch, which is connected to a handle that houses the printed circuit board.

CONSULT THE LOCAL DISTRIBUTOR
OR MANUFACTURER FOR THE AP-
PROPRIATE INTERFACE CABLES.

INDICATIONS FOR USE

The REFSTAR RMT External Reference Patch is indicated for use with the NaviStar RMT Catheter and the CARTO RMT System to provide catheter tip location as well electrogram information.

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CONTRAINDICATIONS

There are no known contraindications for the use of the electrode-patch.

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
INSTRUCTIONS FOR USE

Please refer to these Instructions for Use: the Biosense Webster NAVISTAR RMT Catheter Instructions for Use, and the CARTO RMT System User Manual when using the REFSTAR RMT External Reference-Patch in conjunction with the NAVISTAR RMT Catheter and the CARTO RMT System.

- Remove the external reference patch from its package. Affix the patch to the patient's back (as shown below). Connect the proximal handle of the electrode-patch to the appropriate equipment.
- The REFSTAR RMT External Reference-Patch should be securely positioned on the patient's back. The patch should be placed below the scapula in the area of the seventh (for CARTO RMT System) thoracic vertebra. The body (shaft) of the reference device

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
should follow the long axis of the patient towards the feet or head. Both the patch and the bpd of the reference device should be securely positioned and taped in such a manner as not to subject the device to any abnormal stresses during the procedure which may dislodge the reference device from its location. The handle of the reference device should be connected to the appropriate connection scheme leading to the CARTO RMT System patient interface unit. The CARTO RMT System software program is equipped with a positional reference location feature. This feature is used to ensure that the reference device's tip is within the boundaries of the mapping field.

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PRECAUTIONS

- Do not attempt to use the REFSTAR RMT External Reference Patch prior to completely reading and understanding these Instructions for Use.
- This device is for external use only.
- Store in a cool, dry place. Inspect the packaging and electrode-patch prior to use.
- The REFSTAR RMT External Reference Patch is intended for single patient use only.
- Do not resterilize and reuse.
- This device is packaged for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

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There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the product(s) described herein. Under no circumstances shall Biosense Webster, Inc., or its affiliated companies, be liable for any special, direct, incidental, consequential or other damages other than as expressly provided by specific law.

Without limiting the foregoing, Biosense Webster, Inc. Or its affiliated Companies, shall not be liable for any special, direct, incidental, Consequential, or other damages, arising out of the reuse of any Product(s) labeled for single use or where reuse is prohibited by Applicable law.

Descriptions or specifications appearing in Biosense Webster printed matter, including this publication, are informational only and meant solely to generally describe the product at the time of manufacture and are not made or given as a warranty of the prescribed product in any way.

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- NON-STERILE.
- For single use only.
- Do not use if package is open or damaged.

DESCRIPTION

The REFSTAR or REFSTAR PLUS Electrode-Patch is an integral part of a non-fluoroscopic catheter tip location and electrogram capture technology called CARTO and NOGA. When used with the CARTO and NOGA Systems, the location of the navigation catheter tip is compared to the location of this electrode-patch referred to as a reference patch. This reference device consists of a sensor embedded in an adhesive patch, which is connected to a handle or connector.

CONSULT THE LOCAL DISTRIBUTOR OR MANUFACTURER FOR THE APPROPRIATE INTERFACE OR ADAPTOR CABLES.

INDICATIONS

The REFSTAR or REFSTAR PLUS Electrode-Patch is indicated for use with the Biosense Webster Navigation catheters and the CARTO and NOGA Systems to provide a reference point for catheter tip location.

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CONTRAINDICATIONS

There are no known contraindications for the use of the electrode-patch..

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
INSTRUCTIONS FOR USE

Please refer to these Instructions for Use: the Biosense Webster ESOPHASTAR, NAVISTAR or NOGASTAR catheter's Instructions for Use, and the CARTO or NOGA System's User Manual when using the REFSTAR or REFSTAR PLUS External Reference -Patch in conjunction with the ESOPHASTAR, NAVISTAR or NOGASTAR catheters and the CARTO or NOGA System.


- Remove the external reference patch from its package. Affix the patch to the patient's back (as shown below). Connect the proximal handle or connector of the electrode-patch to the appropriate equipment.

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- The REFSTAR or REFSTAR PLUS External Reference Patch should be securely positioned on the patient's back. The patch should be placed below the scapula in the area of the seventh (for CARTO) or the eighth (for NOGA) thoracic vertebra. The body (shaft) of the reference device should follow the long axis of the patient towards the feet or head. Both the patch and the body of the reference device should be securely positioned and taped in such a manner as not to subject the device to any abnormal stresses during the procedure which may dislodge the reference device from its location. The handle or connector of the reference device should be connected to the appropriate connection scheme leading to the CARTO or NOGA System patient interface unit. The CARTO and NOGA System software program is equipped with a positional reference location feature. This feature is used to ensure that the reference device's tip is within the boundaries of the mapping field.

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WARNINGS

- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.
- Do not expose the external reference patch to organic solvents such as alcohol.
- Do not autoclave the external reference patch.

PRECAUTIONS

- Do not attempt to use the REFSTAR or REFSTAR PLUS External Reference Patch prior to completely reading and understanding these Instructions for Use.
- This device is for external use only.
- Store in a cool, dry, dark place. Inspect the packaging and external reference-patch prior to use.
- The REFSTAR or REFSTAR PLUS Electrode-Patch is intended for single patient use only.
- Do not resterilize and reuse.
- Federal (USA) law restricts this device to sale by or on the order of a physician.
- This device is packaged for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious

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disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

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Disclaimer of Warranty and Limitation of Liability.

There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, ON THE PRODUCT(S)

DESCRIBED HEREIN.. Under no circumstances shall Biosense Webster, Inc. OR ITS AFFILIATED COMPANIES, be liable for any special, direct, incidental, or consequential damages, or other damages other than as expressly provided by specific law.

WITHOUT LIMITING THE FOREGOING, BIOSENSE WEBSTER, INC. OR ITS AFFILIATED COMPANIES, SHALL

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**NOT BE LIABLE FOR ANY SPECIAL,
DIRECT, INCIDENTAL, CONSEQUEN-
TIAL, OR OTHER DAMAGES, ARISING
OUT OF THE REUSE OF ANY PROD-
UCT(S) LABELED FOR SINGLE USE
OR WHERE REUSE IS PROHIBITED
BY APPLICABLE LAW.**

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Medical device is stored at a room temperature, at a place protected from direct sunlight at temperature from +5 до +25 °C and air relative humidity from 20% to 80% (without condensation).

The medical device is transported in cool, dry, dark place, at temperature from

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

State of California)
)
County of Orange)

On November 9, 2017 before me, Richard G. Perez, NOTARY PUBLIC, personally appeared Joanne Whiteman, who proved to me on the basis of satisfactory evidence to be the person whose name is subscribed to the within instrument and acknowledged to me that she executed the same in her authorized capacity and that by her signature on the instrument the person, or the entity upon behalf of which the person acted, executed the instrument.

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WITNESS my hand and official seal.

Signature: 
Richard G. Perez
Notary Public

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Instruction for Use External Reference Patches

- NON-STERILE
- For single use only
- Do not use if package is open or damaged

DESCRIPTION

The CARTO 3 System External Reference Patches allow nonfluoroscopic catheter navigation when used with the CARTO 3 System. The patches provide a means for attachment of the patch sensor cables to the patient in order to convey catheter location. Please refer to the CARTO 3 System's User Manual for additional connectivity information.

INDICATIONS FOR USE

The CARTO 3 System External Reference Patches are indicated for use with the CARTO 3 System to provide a means of attachment of the patch sensor cable to the patient in order to convey catheter location.

CONTRAINDICATIONS

There are no known contraindications for the use of the external reference patches.

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INSTRUCTIONS FOR USE

Please refer to this Instructions for Use and the CARTO 3 System's User Manual when using the CARTO 3 System External Reference Patches in conjunction with the CARTO 3 System.

1. Ensure that the skin surface is completely shaved and dry.
2. Remove each External Reference Patch from the package.
3. Inspect the External Reference Patch to ensure that it is in good physical condition. Care should be taken when opening the packaging to ensure that the sealing mechanism is not damaged.
4. The patches should be placed surrounding the heart and within the CARTO 3 mapping zone. The following is the recommended scheme for patch placement.
5. Back Patch Placement:
 - Attach three (3) patches to the patch sensor cables, and then attach to the patient's back while the patient is sitting.
 - It is recommended to identify the center of the heart, slightly left of the spine, approximately at the level of the 7th thoracic vertebra (T7).
 - Place the center of the first External Reference Patch approximately 12 cm (4.725") superior to the center of the heart.

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- Place the center of the second External Reference Patch approximately 12 cm (4.725") inferior and to the right of the center of the heart.

- Place the center of the third External Reference Patch approximately 12 cm (4.725") inferior and to the left of the center of the heart.

- The back patches should be positioned in a triangle as demonstrated in the Figure 1 below. After the patches have been attached to the patient's back, the patient should lie down.

6. Slide the location pad under the patient's table, such that the three back patches are in the center of the location pad. Verify the position of the patches using the "Location Set Up Phase" on the CARTO 3 System. Adjust the location pad location if required.

7. Chest Patch Placement - The chest patches should be positioned anatomically according to the instructions below. Refer to Figure 2.

- While the patient is lying down, attach three (3) patches to the patch sensor cables.

- Place the first patch on the right lower rib cage (5th / 6th intercostals space) at the right mid-clavicular line. Use the Location Setup phase to verify the patch is within the accurate mapping zone.

- Place second patch on the left lower rib cage (5th / 6th intercostals space) at the anterior axillary line (under V5 ECG patch). Use the Location Setup phase to verify the patch is within the accurate mapping zone.

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- Place the third patch under the clavicle, left or right to the sternum (depending of the Defibrillation Pad position and presence of a device such as ICD or pace maker, etc.). Use the Location Setup phase to verify the patch is within the accurate mapping zone.

8. Using the Location Setup screen, verify that all chest patches are located within the boundaries of the accurate mapping zone.


Figure 1: Back Patch Placement


Figure 2: Chest Patch Placement

For additional education and materials for patch placement, please contact your local Biosense Webster representative.

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WARNINGS

- Discontinue use of this product if skin sensitivity is observed.

PRECAUTIONS

- Do not overlap the External Reference Patches with any other patches on the body including but not limited to EKG patches, defibrillator patches, indifferent electrode patches, etc. or any other External Reference Patch.
- Do not place the External Reference Patches over devices such as pace makers or ICD's.
- Avoid large bony areas such as the scapula and sternum when placing the External Reference Patches.
- Do not expose the external reference patch to organic solvents such as alcohol. If an organic solvent is used on patient, allow to dry for 1 minute prior to applying patches.
- The External Reference Patches should only be connected to the CARTO 3 System.
- Do not attempt to use the External Reference Patches prior to completely reading and understanding these Instructions For Use. This device is for external use only.
- Store in a cool, dry, dark place. Inspect the packaging and external reference patches prior to use.

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- If all External Reference Patches are not used during a procedure, ensure that unused patches are stored in the original, sealed packaging. If an External Reference Patch has not been stored in its original, sealed packaging, the patches should not be used.

- The External Reference Patches are intended for single patient use only.

- Do not sterilize or reuse.

- This device is packaged for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

- If the External Reference Patch does not adhere well to the patient's skin, replace with a new patch.

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WARRANTIES

The manufacturer warrants conformance of this medical device to the declared characteristics within shelf life

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(sterility term) provided that all requirements of the manufacturer's use instructions are met.

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There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the product(s) described herein. Under no circumstances shall Biosense Webster, Inc. or its affiliated companies, be liable for any special, direct, incidental, consequential, or other damages other than as expressly provided by specific law.

WITHOUT LIMITING THE FOREGOING, BIOSENSE WEBSTER, INC. OR ITS AFFILIATED COMPANIES, SHALL NOT BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES, ARISING OUT OF THE REUSE OF ANY PRODUCT(S) LABELED FOR SINGLE USE OR WHERE REUSE IS PROHIBITED BY APPLICABLE LAW.

Descriptions and specifications appearing in Biosense Webster, Inc. printed matter, including this publication, are informational only and meant solely to generally describe the product at the time of manufacture and are not made or given as a warranty of the prescribed product in any way.

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Instruction for Use Catheter-Electrode electrophysiology

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
**I hereby certify accuracy, correctness and reliability
of this document translated into Russian**

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

WITNESS my hand and official seal.

Signature: 
Steve Oberlander
Notary Public

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Instruction for Use Catheter-electrode electrophysiology

- STERILE. Sterilized using ethylene oxide.
- For single use only.
- Do not use if the package is open or damaged.
- Federal (USA) law restricts this device to sale by or on the order of a physician.

DESCRIPTION

The Fixed Curve Catheter-Electrode Electrophysiology (further – catheter) has been designed to facilitate electrophysiological mapping of cardiac structures. The catheter has a high-torque shaft with an array of platinum electrodes at the distal tip that can be used for stimulation and recording. The tip section has a preformed curve with a nylon tip protector. (Standard curve types are identified by a color band on the shaft near the proximal end - see catalog).

The Fixed Curve Lumen Catheter has been designed to facilitate electrophysiological mapping of the heart. The catheter has a high-torque through-lumen shaft with a fixed curve section containing platinum electrodes at the distal tip that can be used for stimulation and recording. The tip section has a preformed curve and is packaged with a nylon tip protector.

The Fixed Curve Lumen Catheter terminates with a luer hub that may be used with a guidewire or for injection purposes.

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The Fixed Orthogonal Catheter has been designed to facilitate electrophysiological mapping of the heart. The catheter has a high-torque shaft with an array of platinum plaque electrodes at the distal tip that can be used for stimulation and recording. The tip section has a preformed curve and is packaged with a nylon tip protector.

The high-torque shaft allows the plane of the distal curve to be manually rotated to assist in positioning the catheter tip at the desired site.

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INDICATIONS FOR USE

The catheter is indicated for electrophysiological mapping of cardiac structures; i.e., stimulation and recording only. In addition, the Fixed Orthogonal Catheter is designed for use in the coronary sinus.

CONTRAINDICATIONS

- The Fixed Curve Catheter has not been shown to be safe and effective for electrical ablation or for use in the coronary vasculature other than the coronary sinus ostium.
- Use of the catheter may not be appropriate for patients with prosthetic valves.

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- A relative contraindication for cardiac catheter procedures is active systemic infection.

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INSTRUCTIONS FOR USE

- Remove the catheter from its package and place it in a sterile work area. Remove nylon tip protector.
- Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
- Connect the lead pins or connectors to the appropriate recording equipment.

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- Advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning. The electrode catheter shaft may be gently hand rotated to facilitate positioning at the desired mapping site.

- Remove the catheter and dispose it in an appropriate manner. Do not resterilize and reuse.

IF THERE ARE ANY QUESTIONS REGARDING THE USE OR PERFORMANCE OF THIS PRODUCT, PLEASE CONSULT WITH THE LOCAL DISTRIBUTOR OR THE MANUFACTURER.

WARNINGS

Cardiac catheterization procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging. Cardiac catheterization should only be performed after adequate attention has been given to potential radiation exposure associated with the procedure, and steps taken to minimize this exposure. Careful consideration must therefore be given for the use of this catheter in pregnant women.

- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.

- Do not expose the catheter to organic solvents such as alcohol.

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PRECAUTIONS

- Do not attempt to operate the Fixed Curve Catheter prior to completely reading and understanding these directions for use.
- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- Store in a cool, dry, dark place.
- The sterile packaging and catheter should be inspected prior to use.
- Catheters are intended for single patient use only.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause

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patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

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There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Product(s) described herein. Under no circumstances shall Biosense Webster, Inc. or its affiliated companies, be liable for any special, direct, incidental, consequential or other damages other than as expressly provided by specific law.

Without limiting the foregoing, Biosense Webster, Inc. or its affiliated companies, shall not be liable for any special, direct, incidental, consequential, or other damages, arising out of the reuse of any

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product(s) labeled for single use or where reuse is prohibited by applicable law.

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Medical device is stored at a room temperature, at a place protected from direct sunlight, at the temperature from +5 to +25 °C and air relative humidity from 20% to 80% (without condensation).

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Medical device in operation is stable against influence of climatic factors for temperature condition from +32 °C to +42 °C and to influence of biological liquids.

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STERILIZATION

- STERILE. Sterilized using ethylene oxide.
- For single use only.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

SHELF LIFE

Use the device prior to the “Use By” date on the package label.

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Instruction for Use Catheters-Electrodes electrophysiology

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
I hereby certify accuracy, correctness and reliability
of this document translated into Russian

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

Instruction for Use Catheters-electrodes electro- physiology

- Sterile. Sterilized using ethylene oxide.
- For single use only.
- Do not use if the package is open or damaged.
- Federal (USA) law restricts this device to sale by or on the order of a physician.

DESCRIPTION

The Catheter-Electrode Electrophysiology Deflectable Tip Electrode (further – catheter) has been designed to facilitate electrophysiological mapping of the heart. The catheter has a high-torque shaft with a deflectable tip section containing an array of platinum electrodes that can be used for stimulation and recording. Several tip configurations are available. (Standard curve types are identified by a color band on the shaft near the proximal end - see catalog.) The Deflectable Tip Orthogonal Catheter has been designed to facilitate electrophysiological mapping of the heart. The catheter has a high-torque shaft with a deflectable tip section containing platinum plaque electrodes that can be used for stimulation and recording.

The Deflectable Braided Tip Catheter is a multi-electrode catheter with a deflectable tip designed to facilitate electrophysiological mapping of the heart. The catheter has a high-torque shaft with a deflectable tip section containing platinum electrodes. The catheter is designed to facilitate stable local electrograms due to the internal

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braid. All of the electrodes may be used for recording and stimulation purposes.

The Deflectable Tip CRISTACATH and ISMUS catheters have been designed for electrophysiological mapping of cardiac structures. The catheters have a high-torque shaft with a tip section containing ten pairs of platinum electrodes that can easily be seen under fluoroscopy. The Deflectable Tip CRISTACATH and ISMUS CATH catheters facilitate simultaneous local electrograms due to the greater number of electrodes and their deflection capabilities.

The model with Auto ID the is equipped with Electronically Erasable Programmable Read Only Memory (EEPROM) which is used to store unique catheter identification information. CARTO® EP Navigation Systems equipped with Auto ID Technology can access the stored information and automatically recognize the catheter information.

The Deflectable Tip HALO XP and HALO catheters have been designed for electrophysiological mapping of the heart, in particular, the tricuspid annulus. The catheters have a high torque shaft with a halo-shaped tip section containing ten pairs of platinum electrodes that can easily be seen under fluoroscopy. The tip section also contains a radiopaque mark in the center of the electrode array. The halo-shaped tip section has a preformed curve which can be positioned around the atrial aspect of the tricuspid annulus. The tip section is packaged with a plastic tip protector.

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The model with AUTO ID is equipped with Electronically Erasable Programable Read Only Memory (EEPROM) which is used to store unique catheter identification information. CARTO® EP Navigation Systems equipped with Auto ID Technology can access the stored information and automatically recognize the catheter information.

A piston in the handpiece is attached to an internal puller wire which changes the radius of curvature. When the piston is pushed forward, the radius of curvature of the preformed loop is reduced; when the thumbknob is pulled back, the radius of the curvature is increased until the tip section returns to the preformed shape. The high torque shaft allows the plane of the curve to be maneuvered in order to facilitate accurate positioning.

The Deflectable Tip HALO XP and HALO catheters facilitate simultaneous local electrograms spanning the tricuspid annulus, from the midseptal to anterior to lateral to posterolateral. Recordings of the entire annulus can be obtained without repositioning the catheter tip.

Tip deflection is controlled at the proximal end by a tubular handpiece in which a piston slides. When the piston is pushed forward with the thumbknob, the tip is deflected (curved). When the piston is pulled back, the tip straightens. The high-torque shaft allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site.

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The catheter interfaces with standard recording equipment via interface cables with the appropriate connectors.

The AUTO ID Models interfaces with CARTO® EP Navigation Systems equipped with Auto ID Technology via interface cables with the appropriate connectors.

CONSULT THE LOCAL DISTRIBUTOR OR MANUFACTURER FOR THE APPROPRIATE INTERFACE CABLES.

INDICATIONS FOR USE

The catheter is indicated for electrophysiological mapping of cardiac structures; i.e., stimulation and recording only. In addition, the Deflectable Tip Orthogonal Catheter is designed for use in the coronary sinus. In addition, the Deflectable Braided Tip Catheter is designed for use in the right side of the heart. In addition, the CRISTACATH and ISMUS catheters are designed to facilitate electrograms in the atrial region of the heart. In addition, the HALO XP and HALO catheters' preformed halo shape of the tip section is designed specifically for the tricuspid annulus.

CONTRAINDICATIONS

The Deflectable Tip Electrode Catheter has not been shown to be safe and effective for electrical ablation or for use in the coronary vasculature other than the coronary sinus

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ostium. In addition, the Braided Tip Catheter is contraindicated for use in the left side of the heart. In addition, the transseptal approach is contraindicated.

Use of the catheter may not be appropriate for patients with prosthetic valves. A relative contraindication for cardiac catheter procedures is active systemic infection.

ADVERSE REACTIONS

A number of serious adverse reactions have been documented for cardiac catheterization procedures including pulmonary embolism, myocardial infarction, stroke, cardiac tamponade, and death. The following complications associated with cardiac catheterization have also been reported in the literature: vascular bleeding, local hematomas, thrombosis, AV fistula, pseudoaneurysm, thromboembolism, vasovagal reactions, cardiac perforation, air embolism, arrhythmias, valvular damage, pneumothorax and hemothorax.

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INSTRUCTIONS FOR USE

- Remove the catheter from its package and place it in a sterile work area. Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.

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- Connect the interface connectors to the appropriate recording equipment or to the Patient Interface Unit of the appropriate CARTO® EP Navigation System and connect the catheter to the interface cable in the models with AUTO ID

- Confirm that the thumbknob is pulled back completely before insertion. Advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning. Adjust the radius of curvature as necessary by manipulating the thumbknob. Pushing the thumbknob forward causes the catheter tip to bend (curve); when the knob is pulled back, the tip straightens.

- Prior to removal of the catheter, confirm that the thumbknob has been pulled back completely. Remove the catheter and dispose of it in an appropriate manner. Do not resterilize and reuse.

- If there are any questions regarding the use or performance of this product, please consult with the local distributor or the manufacturer.

WARNINGS

Cardiac catheterization procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging. Cardiac catheterization should only be performed after adequate attention has been given to potential radiation exposure associated with the procedure, and steps taken to minimize this exposure. Careful consideration

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must therefore be given for the use of this catheter in pregnant women.

- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.
- Do not expose the catheter to organic solvents such as alcohol.
- Do not autoclave the catheter.

PRECAUTIONS

- Do not attempt to operate the Biosense Webster deflectable tip electrophysiology catheter prior to completely reading and understanding these directions for use.
- Store in a cool, dry, dark place.
- The sterile packaging and catheter should be inspected prior to use.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

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- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.

- Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered.

- Always pull the thumbknob of the catheter back before insertion or withdrawal to assure that the catheter tip assumes its original shape.

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There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the product(s) described herein. Under no circumstances shall Biosense Webster, Inc. or its affiliated companies, be liable for

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any special, direct, incidental, consequential or other damages other than as expressly provided by specific law.

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Medical device is stored at a room temperature, at a place protected from direct sunlight at temperature from +5 to + 25 °C and air relative humidity from 20% to 80% (without condensation).

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Medical device in operation is stable against influence of climatic factors for temperature condition from $+32^{\circ}\text{C}$ to $+42^{\circ}\text{C}$ and to influence of biological liquids.

STERILIZATION

- Sterile. Sterilized using ethylene oxide.
- For single use only.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

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Instruction for Use Catheter – Electrode Electrophysiology

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- STERILE. Sterilized using ethylene oxide.
- For single use only.
- Do not use if the package is open or damaged.

DESCRIPTION

The Diagnostic/Ablation electrophysiological Catheter - electrode is a multi-electrode catheter with a deflectable tip designed to facilitate electrophysiological mapping of the heart and to transmit radiofrequency current to the catheter tip electrode for ablation purposes. For ablation, the catheter is used in conjunction with a radiofrequency generator and a dispersive pad (reference electrode). Two models of the catheter are available: the standard model and the temperature sensing model. The temperature sensing model incorporates a temperature sensor in the tip electrode which permits temperature monitoring during radiofrequency ablation. The temperature sensor is embedded in the large tip electrode and electrical connection to the sensor is made via an accessory cable. The catheter has a high-torque shaft with a deflectable tip section containing an array of platinum electrodes. All of the electrodes may be used for recording and stimulation purposes. The large tip electrode also serves to deliver radiofrequency (RF) current from the radiofrequency lesion generator to the desired ablation site. The cylindrical tip electrode may have either straight or concave walls. Tip deflection is controlled at the proximal end by a tubular handpiece in which a piston slides;

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a thumbknob on the piston controls piston travel. When the thumbknob is pushed forward, the tip is deflected (curved). When the thumbknob is pulled back, the tip straightens. The shape of the curve depends on the deflectable tip length (1.5" - 4.0") and the location of the puller wire anchor (proximal or distal) in the deflectable tip. Eleven different curves, designated "A" through "D", posteroseptal D, "E" through "H", "J" and "K" are available. The high-torque shaft also allows the plane of the curved tip to be rotated so that both tip curvature and rotation can assist the physician in accurately positioning the catheter tip in the heart chambers for mapping and radiofrequency ablation. The catheter interfaces with standard recording equipment and a compatible radiofrequency generator via accessory extension cables with the appropriate connectors.

The deflectable braided tip catheter is a temperature sensing catheter that has a high-torque shaft with a deflectable tip section containing platinum electrodes. The catheter is designed to facilitate stable local electrograms due to the internal braid.

The CELSIUS DS Diagnostic/Ablation electrophysiological Deflectable Tip Catheter has an 8 mm deflectable tip with dual temperature sensors designed to facilitate electrophysiological mapping of the heart and to transmit radiofrequency (RF) current to the catheter tip electrode for ablation purposes. For ablation, this catheter is used in conjunction with a compatible RF generator, compatible interface cables and

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a grounding (dispersive) pad. To obtain information about generators that may be compatible with the Biosense Webster Dual Sensor 8 mm catheter, please contact a Biosense Webster representative. Currently, the available curves for these catheters are B, D, F and J; however, future available curves will be all of the standard curves, A-J.

Tip deflection is controlled in the same manner as the standard deflectable tip diagnostic/ablation catheter.

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INDICATIONS FOR USE

The Diagnostic/Ablation electrophysiological Deflectable Tip Catheter is indicated for cardiac electrophysiological mapping (stimulation and recording) and, when used in conjunction with a radiofrequency generator, for cardiac ablation.

THE CELSIUS DS DIAGNOSTIC/ABLATION DEFLECTABLE TIP CATHETER IS INDICATED FOR USE ONLY WITH A COMPATIBLE RF GENERATOR AND APPROPRIATE INTERFACE CABLES. REFER TO THE RF GENERATOR MANUAL FOR PROPER CONNECTION OF THE CATHETER TO THE GENERATOR AND FOR DETAILED INSTRUCTIONS AS TO GENERATOR OPERATION.

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CONTRAINDICATIONS

The catheter should not be used for radiofrequency ablation if the patient has intracardiac mural thrombus or has had a ventriculotomy or atriotomy within the preceding four weeks. The catheter may not be appropriate for patients with prosthetic valves.

A relative contraindication for catheter ablation procedures is active systemic infection.

The transseptal approach is contraindicated in patients with left atrial thrombus or myxoma, or interatrial baffle or patch.

The retrograde trans-aortic approach is contraindicated in patients with aortic valve replacement.

The contraindications for deflectable braided tip catheter are the same as those for the standard deflectable tip diagnostic/ablation catheter. In addition, the 8 French braided tip catheter should not be used in the left side of the heart.

ADVERSE REACTIONS

A number of serious adverse reactions have been documented for catheter ablation procedures including pulmonary embolism,

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myocardial infarction, stroke, cardiac tamponade, and death.

The following complications have been reported in the literature:

- Catheterization/catheter procedure related: vascular bleeding/local hematomas, thrombosis, AV fistula, pseudoaneurysm, thromboembolism and vasovagal reactions, cardiac perforation, tamponade, thrombi, air embolism, arrhythmias and valvular damage, pneumothorax and hemothorax.
- Radiofrequency related: chest pain/discomfort, ventricular tachyarrhythmia, TIA, CVA, complete heart block, coronary artery spasm, coronary artery thrombosis, coronary artery dissection, cardiac thromboembolism, pericarditis, cardiac perforation/tamponade, valvular damage and increased phosphokinase level.

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INSTRUCTION FOR USE

- Remove catheter from package and place in a sterile work area.
- Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
- Connect the catheter to the recording equipment and/or the radiofrequency generator using the appropriate interface cables.

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- Advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning.

The catheter tip can be deflected to facilitate positioning by using the thumbknob to vary tip curvature. Pushing the thumbknob forward causes the catheter tip to bend; when the knob is pulled back, the tip straightens. When it has been determined that the tip electrode is in stable contact with the intended ablation site, the catheter tip electrode connection must be switched from the recording equipment to the radiofrequency generator in preparation for delivery of radiofrequency current. Circuit impedance should be approximately 100 ohms upon initiation of radiofrequency current.

Radiofrequency current may be re-applied to the same or alternate sites using the same catheter. However, in the event of a generator cutoff (impedance or temperature), the catheter must be withdrawn and the tip electrode cleaned of coagulum before radiofrequency current is re-applied. When the catheter is withdrawn, the tip electrode should be cleaned prior to re-insertion. A sterile gauze pad dampened with sterile saline may be used to gently wipe the tip section clean; do not scrub or twist the tip electrode or damage to the tip electrode bond may occur and loosen the tip electrode.

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WARNINGS

- Implantable pacemakers and Automatic Implantable Cardioverter/Defibrillator (AICDs) may be adversely affected by radiofrequency current. It is important to: have temporary external sources of pacing available during ablation; temporarily re-program the pacing system to minimum output or OFF mode to minimize risk of inappropriate pacing. Exercise extreme caution during ablation when in close proximity to atrial or ventricular permanent leads; program the AICD to the OFF mode during the ablation procedure; and, perform complete implantable device analysis on all patients after ablation.
- Patients undergoing septal accessory pathway ablation are at risk for complete AV block which requires the implantation of a permanent pacemaker.
- During the trans-aortic approach, adequate fluoroscopic visualization is necessary to avoid placement of the catheter in the coronary vasculature. Intracoronary placement of the ablation catheter or RF energy application or both have been associated with myocardial infarction.
- Catheter ablation procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging. Catheter ablation should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps

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taken to minimize this exposure. Careful consideration must therefore be given for the use of the device in pregnant women.

- Do not expose catheter to organic solvents such as alcohol.
- Do not autoclave the catheter.
- Do not immerse proximal handle or cable connector in fluids; electrical performance could be affected.
- Do not scrub or twist the distal tip electrode during cleaning.
- Patients undergoing ablations on the septal side of the isthmus may be at risk for complete AV block which requires the implantation of a permanent pacemaker. Permanent pacing may be required in patients who experience inadvertent complete AV block as a result of radiofrequency ablation.
- Do not exceed 70 watts during ablation. If using a 100 watt generator, program the power to cut off at 70 watts.
- This device has not been shown to be safe at temperatures above 60°C.

PRECAUTIONS

- The long-term risks of protracted fluoroscopy and creation of radiofrequency induced lesions have not been established. Careful consideration must therefore be

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given for the use of the device in prepubescent children. Furthermore, the risk/benefit in asymptomatic patients has not been studied.

- To avoid thromboemboli, intravenous heparin is used when entering the left heart during ablation, and many physicians prescribe aspirin, less often warfarin, for about 3 months afterward. No consensus yet exists about the need for short-term anticoagulation after ablation.

- Cardiac ablation procedures should be performed by appropriately trained personnel in a fully-equipped electrophysiology laboratory.

- Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade.

- Catheter advancement should be done under fluoroscopic guidance.

- Do not use excessive force to advance or withdraw the catheter when resistance is encountered.

- Always pull the thumbknob back to straighten the catheter tip before insertion or withdrawal of the catheter.

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- Do not attempt to operate the Diagnostic/Ablation Deflectable Tip Catheter or the radiofrequency generator prior to completely reading and understanding the applicable instructions for use.

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- The Diagnostic/Ablation electrophysiological Deflectable Tip Catheter is intended for use with compatible* radiofrequency generators and accessories only.

- The sterile packaging and catheter should be inspected prior to use.

- Store in a cool, dry, dark place.

- The Diagnostic/Ablation electrophysiology Deflectable Tip Catheter is intended for single patient use only.

- Do not resterilize and reuse.

- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of

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- It is recommended that less than thirty (30), 1 minute applications of energy be performed per catheter.

- The Diagnostic/Ablation Deflectable Tip Catheter used in conjunction with a radiofrequency generator is capable of delivering significant electrical power. Patient or operator injury can result from improper handling of the catheter and dispersive electrode, particularly when operating the device. During energy delivery, the patient should not be allowed to come in contact with grounded metal surfaces. If during ablation, tip temperature does not rise, discontinue ablation immediately and replace catheter.

- The risk of igniting flammable gases or other materials is inherent in electrosurgery. Precautions must be taken to restrict flammable materials from the electrosurgical site.

- Electromagnetic interference (EMI) produced by the Diagnostic/Ablation electrophysiology Deflectable Tip Catheter used in conjunction with a Radiofrequency generator during normal operation may adversely affect the performance of other equipment.

- Electrodes and probes of monitoring and stimulating devices can provide paths for high frequency current. The risk of burns

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can be reduced but not eliminated by placing the electrodes and probes as far away as possible from the ablation site and the dispersive electrode. Protective impedances may reduce the risk of burns, and permit continuous monitoring of the electrocardiogram during energy delivery.

- The temperature sensing version of the catheter measures electrode tip temperature, not tissue temperature. If the generator does not display temperature, verify that the appropriate cable is plugged into the generator. If temperature still is not displayed, there may be a malfunction in the temperature sensing system which must be corrected prior to applying RF power.

- Regularly inspect and test re-usable cables and accessories.

*Laboratory compatibility test on file

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RADIOFREQUENCY ABLATION

For radiofrequency ablation, the catheter must be connected to the appropriate input connector(s) on the radiofrequency generator. For proper RF generator interface, use only a Biosense Webster interface cable.

NOTE: For catheters and cables with tip pin terminations, the pins are marked as follows: a) "D" denotes the tip (ablating)

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electrode, б) “THR” denotes a thermistor temperature sensor, and с) “TC” denotes a thermocouple temperature sensor. The catheter tip electrode must be switched from the electrogram recording equipment (via the generator controls or an external switch) to the radiofrequency generator power output for ablation (see applicable radiofrequency generator manual for instructions). To complete the electrical circuit, a dispersive pad must be connected to the reference electrode input on the generator. Circuit impedance prior to radiofrequency ablation should be approximately 100 ohms. For the temperature sensing catheter, verify that the generator displays a temperature near 37°C after the catheter is inserted in the patient and before applying RF power.

GENERATOR OPERATION

Refer to the applicable radiofrequency generator manual for proper connection of the catheter to the generator and for detailed instructions as to generator operation for radiofrequency ablation.

If radiofrequency current is interrupted for either a temperature or an impedance rise (a preset limit is exceeded), the catheter should be removed, and the tip cleaned of coagulum. When cleaning the tip electrode,

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be careful not to twist the tip electrode with respect to the catheter shaft; twisting may damage the tip electrode bond and loosen the tip electrode.

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There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Biosense Webster, Inc. product(s) described in this publication. Under no circumstances shall Biosense Webster, Inc. or its affiliated companies, be liable for any special, direct, incidental, consequential or other damages other than as expressly provided by specific law.

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Medical device is stored at a room temperature, at a place protected from direct sunlight at the temperature from +5 to +25 °C and air relative humidity from 20% to 80% (without condensation).

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STERILITY

- STERILE. Sterilized using ethylene oxide.
- For single use only.

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Instruction for Use Electrophysiology Catheter- electrode

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Electrophysiology Catheter- electrode

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DESCRIPTION

The EZ STEER DS Bi-Directional Diagnostic/Ablation electrophysiology Deflectable Tip Catheters-electrodes are multi-electrode devices that can be used for electrophysiological mapping and radiofrequency (RF) ablations of the heart. Standard features of this catheter include a braided bi-directional deflectable tip section, an array of four platinum electrodes which includes a 4mm tip dome, a high torque shaft, and a choice of either thermistor or thermocouple temperature sensing capabilities. Additionally, a variety of curve types are available in symmetric or asymmetric combinations, providing two 180° opposed, single planed curves. When used in conjunction with standard recording equipment, a compatible RF generator and a dispersive pad reference electrode all four of the electrodes are capable of intracardiac recording and stimulation while only the tip electrode can deliver RF energy.

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The EZ STEER DS Bi-Directional Diagnostic/Ablation electrophysiological Deflectable Tip Catheter has an 8 mm deflectable tip with dual temperature sensors designed to facilitate electrophysiological mapping of the heart and to transmit radiofrequency (RF) current to the catheter tip electrode for ablation purposes. For ablation, this catheter is used in conjunction with a compatible RF generator, compatible interface cables and a grounding (dispersive) pad.

To obtain information about other generators that may be compatible with the Biosense Webster Dual Sensor 8 mm catheter, please contact a Biosense Webster representative.

Currently, the available curves EZ STEER DS Bi-Directional Catheters include DD, FF, JJ, DF, and FJ. A Rocker Lever is used in the same manner for both models to deflect the tip.

CONSULT THE LOCAL DISTRIBUTOR FOR THE APPROPRIATE INTERFACE CABLES.

INDICATIONS FOR USE

The EZ STEER DS Bi-Directional Diagnostic/Ablation electrophysiological Deflectable Tip Catheters are indicated for cardiac electrophysiological mapping (stimulation and recording) and, when used in conjunction with a radiofrequency generator, for cardiac ablation.

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Detailed Device Description

The EZ STEER DS Bi-Directional model has an 8 mm electrode with three platinum ring electrodes. Each of the ring electrodes is capable of performing both intracardiac recording and stimulation; however, the tip electrode is the only electrode that can be used to deliver RF energy to the desired tissue site. Conductor wires from the electrodes pass through the catheter handle and terminate in the connector.

1. Remove catheter from package and place in a sterile work area.
2. Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
3. Connect the catheter to the recording equipment and/or the radiofrequency generator using the appropriate interface cables.

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
Clockwise rotation from the "off" position increases the friction within the Deflection Mechanism.

Counter-Clockwise rotation from the "on" position decreases the friction within the Deflection Mechanism.

Figure 2

6. When it has been determined that the tip electrode is in stable contact with the intended ablation site, the catheter tip electrode connection must be switched from the recording equipment to the radiofrequency generator in preparation for delivery of radiofrequency current. Circuit impedance should be approximately 100 ohms upon initiation of radiofrequency current.

7. Radiofrequency current may be re-applied to the same or alternate sites using the same catheter. However, in the event of a generator cutoff (impedance or temperature), the catheter must be withdrawn and the tip electrode cleaned of coagulum before radiofrequency current is re-applied. When the catheter is withdrawn, the tip electrode should be cleaned prior to re-insertion. A sterile gauze pad dampened with sterile saline may be used to gently wipe the tip section clean; do not scrub or twist the tip electrode or damage to the tip electrode bond may occur and loosen the tip electrode.

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CURVE TYPES

Currently, the available curves for EZ STEER DS Bi-Directional catheters include DD, FF, JJ, DF, and FJ (Figure 3).

Figure 3

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The EZ STEER DS Bi-Directional Diagnostic/Ablation electrophysiological Deflectable Tip Catheter should be used only with an RF generator that has been shown to be safe and effective for cardiac ablation. Refer to the RF generator manual for detailed generator operating instructions for RF catheter ablation. If radiofrequency current is interrupted for either a temperature or an impedance rise (a preset limit is exceeded), the catheter should be removed, and the tip cleaned of coagulum. When cleaning the tip electrode, be careful not to twist the tip electrode with respect to the catheter shaft; twisting may damage the tip electrode bond and loosen the tip electrode. Use the appropriate Biosense Webster accessory cable to connect the EZ STEER DS Bi-Directional Diagnostic/Ablation electrophysiological Deflectable Tip Catheter to a compatible RF generator.

If there are any questions regarding the use or performance of this product, please consult with the local distributor or the manufacturer.

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WARNINGS

- Significant x-ray exposure, can result in acute radiation injury as well as dose-related risk for somatic and genetic effects (dose = duration of the fluoroscopic imaging multiplied by x-ray beam intensity). Take all appropriate measures to minimize x-ray exposure to both patients and clinical staff.
- Pregnancy -- Careful consideration should be given to the use of this device in pregnant women.
- Ablation from within a coronary artery can cause myocardial injury and death. Adequate fluoroscopic visualization is necessary during the transaortic approach to avoid placement of the ablation catheter in the coronary vasculature.
- Stroke or myocardial infarction may occur in patients undergoing left-sided ablation procedures. Patients should be closely monitored during the post-ablation period for clinical manifestations of embolic events.
- Implantable pacemakers and implantable cardioverter/defibrillators (ICDs) may be adversely affected by RF ablation. ICDs should be deactivated during ablation. Have temporary external sources of pacing and defibrillation available during ablation. Exercise extreme caution during ablation when in close proximity to device leads and perform a complete analysis of the implanted device function after ablation.

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- Complete AV block can occur when ablating septal accessory pathways or in the treatment of AVNRT. Closely monitor AV conduction during RF energy delivery and immediately terminate energy delivery if partial or complete AV block is observed. Using catheters with distal pair electrode spacing greater than 2 mm may increase the risk of AV nodal damage.
- To reduce the risk of brachial plexus injury, physicians are advised to position the patient's arms inferiorly in the normal position with the hands down by the hips as opposed to superiorly.
- Always place the Rocker Lever in the neutral position to straighten the catheter tip before insertion or withdrawal of the catheter.
- Patients undergoing ablations on the septal side of the isthmus may be at risk for complete AV block, which requires the implantation of a permanent pacemaker. Permanent pacing may be required in patients who experience inadvertent complete AV block as a result of radiofrequency ablation.
- Do not exceed 70 watts during ablation. If using a 100-watt generator, program the power to cut off at 70 watts.

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- This device has not been shown to be safe at temperatures above 60° C.

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PRECAUTIONS

- Cardiac ablation procedures should be performed only by appropriately trained personnel in a fully-equipped electrophysiology laboratory.
- Do not attempt to operate the EZ STEER DS Bi-Directional Catheter or the RF generator prior to completely reading and understanding the applicable directions for use.
- The long-term risks of protracted fluoroscopy have not been established. Careful consideration must therefore be given for the use of the device in prepubescent children.
- The long-term risks of lesions created by RF ablation have not been established. In particular, any long-term effects of lesions in proximity to the specialized conduction system or coronary vasculature are unknown.
- It is recommended that no more than seventy (70), one minute applications of energy be performed per catheter.

Catheter Compatibility

- The EZ STEER DS Bi-Directional Catheters are intended for use with a compatible RF generator (see attached supplement pages) and Biosense Webster accessories only.
- Read and follow the dispersive electrode manufacturer's instructions for use; the use of dispersive electrodes,

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which meet or exceed ANSI/AAMI requirements (HF18), is recommended.

Handling and Sterilization

- FOR SINGLE USE ONLY.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Observe "Use By" Date.
- Sterilized with ethylene oxide gas.
- The sterile packaging and catheter should be inspected prior to use. If the package or the catheter appears damaged, do not use. Contact your local Biosense Webster representative.
- Catheter damage may occur due to:
 - autoclaving
 - resterilizing
 - exposure to organic solvents
 - immersing proximal handle or cable connector in fluids

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- Catheter materials are not compatible with magnetic resonance imaging (MRI).
- During ablation procedures, this catheter is used in conjunction with an RF generator. Electromagnetic interference (EMI) produced by the radiofrequency (RF) generator during the delivery of RF power may adversely affect the performance of other equipment.

Precautions During Catheter Use

- The patient should not contact grounded metal surfaces. Use only isolated amplifiers, pacing equipment, and ECG equipment or patient injury death may occur. Leakage current from any connected device to the patient must not exceed 10 microAmps (μA) under any circumstances.

- Do not use excessive force to advance or withdraw the catheter. Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade.
- Do not insert or withdraw the catheter without placing the catheter tip in neutral position (midway between the two curves).
- Use both fluoroscopy and electrograms to monitor the advancement of the catheter to the area of the endocardium under investigation to avoid vascular or cardiac damage.

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RADIOFREQUENCY ABLATION

For radiofrequency ablation, the catheter must be connected to the appropriate input connector(s) on the radiofrequency (RF) generator. For proper RF generator interface, use only a Biosense Webster interface cable. The catheter tip electrode must be switched from the electrogram recording equipment (via the generator controls or an external switch) to the radiofrequency generator power output for ablation (see applicable radiofrequency generator manual for instructions). To complete the electrical circuit, a dispersive pad must be connected to the reference electrode input on the generator. Circuit impedance prior to radiofrequency ablation should be approximately 100 ohms. For the temperature sensing catheter, verify that the generator displays a temperature near 37° C after the catheter is inserted into the patient and before applying RF power.

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GENERATOR OPERATION

Refer to the applicable radiofrequency generator manual for proper connection of the catheter to the generator and for detailed instructions as to generator operation for radiofrequency ablation. Based on historical data obtained from prior clinical studies on similar design catheters, the typical radiofrequency ablation application parameters are 30 watts for 35 seconds when using the output of a 500 kHz radiofrequency generator (Radionics RFG-3D) into a 102 ohm circuit impedance. These parameters are considered typical and will vary depending on the ablation site, the specific conditions present in each procedure and the radiofrequency generator control circuitry. If temperature control is being used, only the desired temperature need be set; the generator will control the power to maintain the set temperature. If radiofrequency current is interrupted for either a temperature or an impedance rise (a preset limit is exceeded), the catheter should be removed, and the tip cleaned of coagulum. When cleaning the tip electrode, be careful not to twist the tip electrode with respect to the catheter shaft; twisting may damage the tip electrode bond and loosen the tip electrode.

WARRANTIES

The manufacturer warrants conformance of this medical device to the declared charac-

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teristics within shelf life (sterility term) provided that all requirements of the manufacturer's use instructions are met.

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There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Biosense Webster, Inc. product(s) described in this publication. Under no circumstances shall Biosense Webster, Inc. or its affiliated companies, be liable for any special, direct, incidental, consequential or other damages other than as expressly provided by specific law. Without limiting the foregoing, Biosense Webster, Inc. or its affiliated companies, shall not be liable for any special, direct, incidental, consequential, or other damages, arising out of the reuse of any product(s) labeled for single use or where reuse is prohibited by applicable law.

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CONDITIONS OF TRANSPORTA- TION, STORAGE AND OPERATION

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This device is packaged and sterilized for single use only. Do not reuse, reprocess or re-sterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or crossinfection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient

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Instruction for Use Catheter-Electrode Electrophysiology

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DESCRIPTION

The WEBSTER Fixed Curve Electrophysiology (EP) Catheter-Electrode (further – Catheter) has been designed to facilitate electrophysiological mapping of cardiac structures. The catheter has a high-torque shaft with an array of platinum electrodes at the distal tip that can be used for stimulation and recording. The tip section has a preformed curve with a nylon tip protector. (Standard curve types are identified by a color band on the shaft near the proximal end - see catalog). The basic catheter interfaces with standard recording equipment via interface cables with the appropriate connects. The catheter with AUTO ID Technology is equipped with Electronically Erasable Programmable Read Only Memory (EEPROM) which is used to store unique catheter identification information. CARTO EP Navigation Systems equipped with Auto ID Technology can access the stored information and automatically recognize the catheter information. The catheter interfaces with CARTO EP Navigation Systems equipped with Auto ID Technology via interface cables with the appropriate connectors.

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CONSULT THE LOCAL DISTRIBUTOR OR MANUFACTURER FOR THE APPROPRIATE INTERFACE CABLES. REFER TO THE APPROPRIATE CARTO EP NAVIGATION SYSTEM USER MANUAL FOR DETAILS ON SYSTEM SETTINGS RELATED TO THE CATHETER.

INDICATIONS FOR USE

The catheter is indicated for electrophysiological mapping of cardiac structures; i.e., stimulation and recording only.

CONTRAINDICATIONS

- The WEBSTER Fixed Curve EP Catheter has not been shown to be safe and effective for electrical ablation or for use in the coronary vasculature other than the coronary sinus ostium.
- Use of the catheter may not be appropriate for patients with prosthetic valves.
- A relative contraindication for cardiac catheter procedures is active systemic infection.

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INSTRUCTIONS FOR USE

- Remove the catheter from its package and place it in a sterile work area. Remove nylon tip protector.
- Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
- Connect the lead pins or connectors to the appropriate recording equipment, or in case of the model with AUTO ID, connect the interface cable to the Patient Interface Unit of

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the appropriate CARTO EP Navigation System and connect the catheter to the interface cable.

- Advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning. The electrode catheter shaft may be gently hand rotated to facilitate positioning at the desired mapping site.
- Remove the catheter and dispose it in an appropriate manner. Do not resterilize and reuse.
- If there are any questions regarding the use or performance of this product, please consult with the local distributor or the manufacturer

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- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.
- Do not expose the catheter to organic solvents such as alcohol.

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PRECAUTIONS

- Do not attempt to operate the WEBSTER Fixed Curve EP Catheter prior to completely reading and understanding these directions for use.
- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- Store in a cool, dry place.
- The sterile packaging and catheter should be inspected prior to use.
- WEBSTER Fixed Curve EP Catheter is intended for single patient use only.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or

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resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

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Without limiting the foregoing, Biosense Webster, Inc. or its affiliated companies,

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shall not be liable for any special, direct, incidental, consequential, or other damages, arising out of the reuse of any product(s) labeled for single use or where reuse is prohibited by applicable law.

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Medical device in operation are stable

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against influence of climatic factors for temperature condition from +32 °C to +42 °C and to influence of biological liquids.

STERILIZATION

- STERILE. Sterilized with ethylene oxide gas.
- For single use only.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

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Instruction for Use Catheter- Electrophysiology

- Sterile. Sterilized with ethylene oxide gas.
- For single use only.
- Do not use if the package is opened or damaged.
- Federal Law (U.S.A) restricts the sale of this device by or on the order of a physician only.

DESCRIPTION

The WEBSTER Duo-Decapolar Catheter-Electrode Electrophysiology (further - catheter) has been designed to facilitate electrophysiological mapping of the heart. The catheter has a high-torque shaft with a deflectable tip section containing an array of platinum electrodes that can be used for stimulation and recording.

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INDICATIONS FOR USE

The catheter is indicated for electrophysiological mapping of cardiac structures; i.e., stimulation and recording only. In addition, the WEBSTER Duo-Decapolar Catheter is designed to facilitate electrogram mapping in the atrial region of the heart and the coronary sinus.

CONTRAINDICATIONS

The WEBSTER Duo-Decapolar Catheter has not been shown to be safe and effective for electrical ablation. In addition, the transseptal approach is contraindicated.

- Use of the catheter may not be appropriate for patients with prosthetic valves.
- A relative contraindication for cardiac catheter procedures is active systemic infection.
- Do not use this device via the retrograde approach in patients with aortic valve replacement.
- Use of the catheter is contraindicated in patients with totally obstructed Coronary Sinus.
- The catheter should not be used in patients unable to receive heparin or an acceptable alternative to achieve the adequate anticoagulation.
- Electrophysiology studies are contraindicated when reversible factors make the

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findings unrepresentative of the patient's disease state (e.g. electrolyte imbalance).

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INSTRUCTIONS FOR USE

- Remove the catheter from its package and place it in a sterile work area. Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
- Connect the interface connectors to the appropriate recording equipment.
- Confirm that the thumbknob is pulled back completely before insertion. Advance the catheter to the area of the endocardium under

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investigation. Use both fluoroscopy and electrograms to aid in proper positioning. Adjust the radius of curvature as necessary by manipulating the thumbknob. Pushing the thumbknob forward causes the catheter tip to bend (curve); when the knob is pulled back, the tip straightens.

- Prior to removal of the catheter, confirm that the thumbknob has been pulled back completely. Remove the catheter and dispose of it in an appropriate manner. Do not resterilize and reuse.
- If there are any questions regarding the use or performance of this product, please consult with the local distributor or the manufacturer.

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- Do not autoclave the catheter.

PRECAUTIONS

- Do not attempt to operate the WEBSTER Duo-Decapolar Catheter prior to completely reading and understanding these directions for use.
- Store in a cool, dry place.
- The sterile packaging and catheter should be inspected prior to use.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.

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incidental, consequential, or other damages, arising out of the reuse of any product(s) labeled for single use or where reuse is prohibited by applicable law.

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STERILIZATION

- Sterile. Sterilized with ethylene oxide gas.
- For single use only.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another.

Contamination of the device may lead to injury, illness, or death of the patient.

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Instruction for Use Catheters-electrodes diagnostic

- Sterile. Sterilized using ethylene oxide.
- For single use only.
- Do not use if the package is open or damaged.
- Federal law (U.S.A) restricts this device to sale by or on the order of a physician.

DESCRIPTION

The Catheter-Electrode Diagnostic Deflectable Tip Electrode (further – catheter) has been designed to facilitate electrophysiological mapping of the heart. The catheter has a high-torque shaft with a deflectable tip section containing an array of platinum electrodes that can be used for stimulation and recording. Several tip configurations are available. (Standard curve types are identified by a color band on the shaft near the proximal end - see catalog.) The Deflectable Tip Orthogonal Catheter has been designed to facilitate electrophysiological mapping of the heart. The catheter has a high-torque shaft with a deflectable tip section containing platinum plaque electrodes that can be used for stimulation and recording.

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INDICATIONS FOR USE

The catheter is indicated for electrophysiological mapping of cardiac structures; i.e., stimulation and recording only. In addition, the Deflectable Tip Orthogonal Catheter is designed for use in the coronary sinus. In addition the Deflectable Braided Tip Catheter is designed for use in the right side of the heart.

CONTRAINDICATIONS

The Deflectable Tip Electrode Catheter has not been shown to be safe and effective for electrical ablation or for use in the coronary vasculature other than the coronary sinus ostium. In addition, the Braided Tip Catheter is contraindicated for use in the left side of the heart. In addition, the transseptal approach is contraindicated.

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INSTRUCTIONS FOR USE

- Remove the catheter from its package and place it in a sterile work area. Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
- Connect the interface connectors to the appropriate recording equipment.
- Confirm that the thumbknob is pulled back completely before insertion. Advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and

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electrograms to aid in proper positioning. Adjust the radius of curvature as necessary by manipulating the thumbknob. Pushing the thumbknob forward causes the catheter tip to bend (curve); when the knob is pulled back, the tip straightens.

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- Do not autoclave the catheter.

PRECAUTIONS

- Do not attempt to operate the deflectable tip electrophysiology catheter prior to completely reading and understanding these directions for use.
- Store in a cool, dry, dark place.
- The sterile packaging and catheter should be inspected prior to use.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter

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advancement and placement should be done under fluoroscopic guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered.

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Without limiting the foregoing, Biosense Webster, Inc. or its affiliated companies, shall not be liable for any special, direct, incidental, consequential, or other damages, arising out of the reuse of any product(s) labeled for single use or

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where reuse is prohibited by applicable law.

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A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy or validity of that document.


of this document translated into Russian

State of California)
)
County of Orange)


On October 19, 2017 before me, Steve Oberlander, NOTARY PUBLIC, personally appeared Joanne Whiteman, who proved to me on the basis of satisfactory evidence to be the person whose name is subscribed to the within instrument and acknowledged to me that she executed the same in her authorized capacity and that by her signature on the instrument the person, or the entity upon behalf of which the person acted, executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

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Instruction for Use Catheter-Electrode Diagnostic

- STERILE. Sterilized with ethylene oxide gas.
- For single use only.
- Do not use if the package is open or damaged.

DESCRIPTION

The NAVISTAR Catheter-Electrode diagnostic (further – catheter) is a multi-electrode Catheter with a deflectable tip designed to facilitate electrophysiological mapping of the heart and to transmit radiofrequency current to the Catheter tip electrode for ablation purposes. For ablation, the Catheter is used in conjunction with a radiofrequency generator and a dispersive pad (reference electrode).

The NAVISTAR DS Diagnostic/Ablation Deflectable Tip Catheter has an 8 mm deflectable tip with dual temperature sensors, and is designed to facilitate electrophysiological mapping of the heart and to transmit radiofrequency (RF) current to the Catheter tip electrode for ablation purposes. For ablation, the Catheter is used in conjunction with a compatible RF generator and a dispersive pad (reference electrode). Currently, the available curves for these Catheters are B, C, D, E, and F; however, future available curves will be all of the standard curves, A-F.

Tip deflection is controlled at the proximal end by a tubular handpiece in which a piston slides. When the piston is pushed forward with the thumbknob, the tip is deflected (curved). When the piston is pulled back, the tip straightens. The high-torque shaft allows the plane of the curved tip to

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be rotated to facilitate accurate positioning of the Catheter tip at the desired site. The Catheter interfaces with standard recording equipment and a compatible radiofrequency generator via accessory extension cables with the appropriate connectors.

The NAVISTAR and NAVISTAR DS Catheters are also an integral part of a non-fluoroscopic Catheter tip location and electrogram capture technology (CARTO EP Navigation System).

For use in non-fluoroscopic mapping procedures, refer to the instructions for the CARTO System. An appropriate location reference device is required for location reference position purposes.

For use in radiofrequency ablation procedures, refer to the radiofrequency generator operating instructions.

INDICATIONS FOR USE

The NAVISTAR and NAVISTAR DS Catheters are indicated for cardiac electrophysiological mapping (stimulation and recording) and, when used in conjunction with a radiofrequency generator, for cardiac ablation.

THE NAVISTAR DS DIAGNOSTIC/ABLATION DEFLECTABLE TIP CATHETER IS INDICATED FOR USE ONLY WITH A COMPATIBLE RF GENERATOR AND APPROPRIATE INTERFACE CABLES. REFER TO THE RF GENERATOR MANUAL FOR PROPER CONNECTION OF THE CATHETER TO THE GENERATOR

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AND FOR DETAILED INSTRUCTIONS AS TO GENERATOR OPERATION OF RF ABLATION.

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CONTRAINDICATIONS

The Catheter should not be used for radiofrequency ablation if the patient has intracardiac mural thrombus or has had a ventriculotomy or atriotomy within the preceding four weeks. The Catheter may not be appropriate for patients with prosthetic valves.

A relative contraindication for Catheter ablation procedures is active systemic infection.

The transseptal approach is contraindicated in patients with left atrial thrombus or myxoma, or interatrial baffle or patch. The retrograde transaortic approach is contraindicated in patients with aortic valve replacement.

ADVERSE REACTIONS

A number of serious adverse reactions have been documented for Catheter ablation procedures including pulmonary embolism, myocardial infarction, stroke, cardiac tamponade, and death.

The following complications were also noted to have occurred during the study or have been reported in the literature:

- Catheterization/Catheter procedure related: vascular bleeding/local hematomas, thrombosis, AV fistula, pseudoaneurysm, thromboembolism and vasovagal reactions, cardiac perforation, tamponade, thrombi, air embolism, arrhythmias and

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valvular damage, pneumothorax and hemothorax.

• Radiofrequency related: chest pain/discomfort, ventricular tachyarrhythmia, TIA, CVA, complete heart block, coronary artery spasm, coronary artery thrombosis, coronary artery dissection, cardiac thromboembolism, pericarditis, cardiac perforation/tamponade, valvular damage and increased phosphokinase level.

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CONDITIONS OF USE

This medical device is used only in the conditions of medical and therapeutic institutions by qualified medical personnel.

INSTRUCTIONS FOR USE

1. Remove Catheter from package and place in a sterile work area.
2. Create a vascular access in a large central vessel using aseptic techniques and insert the Catheter.
3. Connect the Catheter to the recording equipment and/or the radiofrequency generator using the appropriate interface cables.
4. Advance the Catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning.
5. The Catheter tip can be deflected to facilitate positioning by using the thumbknob to vary tip curvature. Pushing the thumbknob forward causes the Catheter tip to bend; when the knob is pulled back, the tip straightens.

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6. When it has been determined that the tip electrode is in stable contact with the intended ablation site, the Catheter tip electrode connection must be switched from the recording equipment to the radiofrequency generator in preparation for delivery of radiofrequency current. Circuit impedance should be approximately 100 ohms upon initiation of radiofrequency current.

7. Radiofrequency current may be re-applied to the same or alternate sites using the same Catheter. However, in the event of a generator cutoff (impedance or temperature), the Catheter must be withdrawn and the tip electrode cleaned of coagulum before radiofrequency current is re-applied. When the Catheter is withdrawn, the tip electrode should be cleaned prior to re-insertion. A sterile gauze pad dampened with sterile saline may be used to gently wipe the tip section clean; do not scrub or twist the tip electrode or damage to the tip electrode bond may occur and loosen the tip electrode.

WARNINGS

- Significant x-ray exposure, can result in acute radiation injury as well as dose-related risk for somatic and genetic effects (dose = duration of the fluoroscopic imaging X x-ray beam intensity). Take all appropriate measures to minimize x-ray exposure to both patients and clinical staff.

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- Pregnancy – Careful consideration should be given to the use of this device in pregnant women.
- Ablation from within a coronary artery can cause myocardial injury and death. Adequate fluoroscopic visualization is necessary during the transaortic approach to avoid placement of the ablation Catheter in the coronary vasculature.
- Stroke or myocardial infarction may occur in patients undergoing left-sided ablation procedures. Patients should be closely monitored during the post-ablation period for clinical manifestations of embolic events.
- Implantable pacemakers and implantable cardioverter/defibrillators (ICDs) may be adversely affected by RF ablation. ICDs should be deactivated during ablation. Have temporary external sources of pacing and defibrillation available during ablation. Exercise extreme caution during ablation when in close proximity to device leads and perform a complete analysis of the implanted device function after ablation.
- Complete AV block can occur when ablating septal accessory pathways or in the treatment of AVNRT. Closely monitor AV conduction during RF energy delivery and immediately terminate energy delivery if partial or complete AV block is observed. Using Catheters with distal pair electrode spacing greater than 2 mm may increase the risk of AV nodal damage.
- To reduce the risk of brachial plexus injury, physicians are advised to position the

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patient's arms inferiorly in the normal position with the hands down by the hips as opposed to superiorly.

The warnings for NAVISTAR DS Catheter are the same as the standard NAVISTAR Diagnostic/Ablation Deflectable Tip Catheter with the additional warnings below.

- Patients undergoing ablations on the septal side of the isthmus may be at risk for complete AV block which requires the implantation of a permanent pacemaker. Permanent pacing may be required in patients who experience inadvertent complete AV block as a result of RF ablation.
- Do not exceed 70 watts during ablation. If using a 100-watt generator, program the power to cut off at 70 watts.
- This device has not been shown to be safe at temperatures above 60°C.

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RADIOFREQUENCY ABLATION

For radiofrequency ablation, the Catheter must be connected to the appropriate input connector(s) on the radiofrequency (RF) generator. For proper RF generator interface, use only a Biosense Webster interface cable.

The Catheter tip electrode must be switched from the electrogram recording

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equipment (via the generator controls or an external switch) to the radiofrequency generator power output for ablation (see applicable radiofrequency generator manual for instructions). To complete the electrical circuit, a dispersive pad must be connected to the reference electrode input on the generator. Circuit impedance prior to radiofrequency ablation should be approximately 100 ohms. For the temperature sensing Catheter, verify that the generator displays a temperature near 37° C after the Catheter is inserted into the patient and before applying RF power.

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GENERATOR OPERATION

Refer to the applicable radiofrequency generator manual for proper connection of the Catheter to the generator and for detailed instructions

as to generator operation for radiofrequency ablation.

Based on historical data obtained from prior clinical studies on similar design Catheters, the typical radiofrequency ablation application parameters are 30 watts for 35 seconds when using the output of a 500kHz radiofrequency generator (Radionics RFG-3D) into a 102 ohm circuit impedance. These parameters are considered typical and will vary depending on the ablation site, the specific conditions present in each procedure and the radiofrequency generator control circuitry. If temperature control is being used, only the desired temperature need be set; the generator will control the power to maintain the set temperature.

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If radiofrequency current is interrupted for either a temperature or an impedance rise (a pre-set limit is exceeded), the Catheter should be removed, and the tip cleaned of coagulum. When cleaning the tip electrode, be careful not to twist the tip electrode with respect to the Catheter shaft; twisting may damage the tip electrode bond and loosen the tip electrode.

PRECAUTIONS

- Cardiac ablation procedures should be performed only by appropriately trained personnel in a fully-equipped electrophysiology laboratory.
- Do not attempt to operate the NAVISTAR or NAVISTAR DS Catheter or the RF generator prior to completely reading and understanding the applicable directions for use.
- The long-term risks of protracted fluoroscopy have not been established. Careful consideration must therefore be given for the use of the device in prepubescent children.
- The long-term risks of lesions created by RF ablation have not been established. In particular, any long-term effects of lesions in proximity to the specialized conduction system or coronary vasculature are unknown.
- When using the NAVISTAR or NAVISTAR DS Catheters with conventional systems (using fluoroscopy to determine Catheter tip location), or with the CARTO System, careful Catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Do not use excessive force to advance or withdraw the Catheter when resistance is encountered.

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- It is recommended that less than thirty (30), 1 minute applications of energy be performed per Catheter.

- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

Catheter Compatibility

- The NAVISTAR and NAVISTAR DS Catheters are intended for use with a compatible RF generator (see attached supplement pages) and Biosense Webster accessories only.

- Read and follow the dispersive electrode manufacturer's instructions for use; the use of dispersive electrodes, which meet or exceed ANSI/AAMI requirements (HF18), is recommended.

Handling and Sterilization

- SINGLE USE ONLY.
- Store in a cool, dry place.
- Observe "Use By" Date.
- Sterilized with ethylene oxide gas.

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- The sterile packaging and Catheter should be inspected prior to use. If the package or the Catheter appears damaged, do not use. Contact your local Biosense Webster representative.

- Catheter damage may occur due to:

- autoclaving
- re-sterilizing
- exposure to organic solvents

- immersing proximal handle or cable connector in fluids

Environmental and EMI

- Catheter materials are not compatible with magnetic resonance imaging (MRI).
- During ablation procedures, this Catheter is used in conjunction with an RF generator. Electromagnetic interference (EMI) produced by the radiofrequency (RF) generator during the delivery of RF power may adversely affect the performance of other equipment.

Precautions During Catheter Use

- The patient should not contact grounded metal surfaces. Use only isolated amplifiers, pacing equipment, and ECG equipment or patient injury or death may occur. Leakage current from any connected device to the patient must not exceed 10 microAmps (μ A) under any circumstances.
- Do not use excessive force to advance or withdraw the Catheter. Careful Catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade.
- Always pull the thumbknob back to straighten the Catheter tip before insertion or withdrawal of the Catheter.

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- Use both fluoroscopy and electrograms to monitor the advancement the Catheter to the area of the endocardium under investigation to avoid vascular or cardiac damage.

Precautions During Ablation

- Do not increase power before checking for lead connection and appropriate dispersive electrode application. Effective contact between the patient and the dispersive electrode must be verified whenever the patient is repositioned.
- Do not deliver RF energy with Catheter outside the target site. The RF generator can deliver significant electrical energy and may cause patient or operator injury.
- Avoid use of electrodes and probes of monitoring and stimulating devices which could provide paths for high frequency current. Reduce the burn hazard by placing the electrodes and probes as far away as possible from the ablation site and the dispersive electrode.
- In the event of a generator cutoff (impedance or temperature), the Catheter must be withdrawn and the tip electrode cleaned of coagulum before RF current is re-applied. Use only sterile saline and gauze pad to clean the tip.
- Do not scrub or twist the tip electrode as damage may cause Catheter failure or patient injury.
- Discontinue ablation immediately and replace Catheter if tip temperature fails to rise during ablation.

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- The Catheter measures electrode tip temperature, not tissue temperature. If the generator does not display temperature, verify that the appropriate cable is plugged into the generator. If temperature still is not displayed, there may be a malfunction in the temperature sensing system which must be corrected prior to applying RF power.

WARRANTIES

The manufacturer warrants conformance of this medical device to the declared characteristics within shelf life (sterility term) provided that all requirements of the manufacturer's use instructions are met.

Disclaimer of Warranty and Limitation of Liability

There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the product(s) described herein. Under no circumstances shall Biosense Webster, Inc. or its affiliated companies, be liable for any special, direct, incidental, consequential or other damages other than as expressly provided by specific law.

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Without limiting the foregoing, Biosense Webster, Inc. or its affiliated companies, shall not be liable for any special, direct, incidental, consequential, or other damages, arising out of the reuse of any product(s) labeled for single use or where reuse is prohibited by applicable law.

Descriptions or specifications appearing in Biosense Webster, Inc. printed matter, including this publication, are informational only and meant solely to generally describe the product at the time of manufacture and are not made or given as a warranty of the prescribed product in any way.

CONDITIONS OF TRANSPORTA- TION, STORAGE AND OPERA- TION

The medical device is transported using all available modes of transportation in the covered vehicles according to the transportation rules forced on each transportation mode.

Medical device is stored at a room temperature, at a place protected from direct sunlight at the temperature from +5 to + 25 °C and air relative humidity from 20% to 80% (without condensation).

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The medical device is transported in cool, dry, dark place, at temperature from -40°C to $+54^{\circ}\text{C}$ and air relative humidity from 20% to 80% (without condensation).

Medical device in operation are stable against influence of climatic factors for temperature condition from $+32^{\circ}\text{C}$ to $+42^{\circ}\text{C}$ and to influence of biological liquids.

STERILIZATION

- STERILE. Sterilized with ethylene oxide gas.
- For single use only.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

SHELF LIFE

Use the device prior to the "Use By" date on the package label.

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Instruction for Use Catheter-Electrode Diagnostic Deflectable tip Catheter

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DESCRIPTION

The NOGASTAR deflectable tip electrode catheter-electrode diagnostic (further — catheter) has been designed to facilitate electromechanical mapping of the heart. The catheter has a high-torque shaft with a deflectable tip section containing an array of platinum electrodes that can be used for recording.

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The NOGASTAR Catheter is an integral part of a nonfluoroscopic catheter tip location and electrogram capture technology known as the NOGA System.

For use in non-fluoroscopic mapping procedures, refer to the instructions for the NOGA System.

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INDICATIONS FOR USE

The NOGASTAR Catheter is indicated for electromechanical mapping of cardiac structures; i.e., recording. In addition, the NOGASTAR Catheter is used with the NOGA System to provide catheter tip location as well as electro-anatomical information.

CONTRAINDICATIONS

Use of the catheter may not be appropriate for patients with prosthetic valves. Do not use this device in patients with active systemic infection.

The transseptal approach is contraindicated in patients with left atrial thrombus or myxoma, or interatrial baffle or patch. Further, the catheter is non-over-the-wire and does not utilize a transseptal sheath system.

ADVERSE REACTIONS

A number of serious adverse reactions have been documented for cardiac catheterization procedures including pulmonary embolism, myocardial infarction, stroke,

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cardiac tamponade, and death. The following complications associated with cardiac catheterization have also been reported in the literature: vascular bleeding, local hematomas, thrombosis, AV fistula, pseudoaneurysm, thromboembolism, vasovagal reactions, cardiac perforation, air embolism, arrhythmias, valvular damage, pneumothorax and hemothorax.

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INSTRUCTIONS FOR USE

(Please refer to both these Instructions for Use and the NOGA System's User Manual when using the NOGASTAR Catheter in conjunction with the NOGA System.)

- Remove the catheter from its package and place it in a sterile work area. Connect the interface connectors to the appropriate recording equipment.
- Create a vascular access in a large central vessel using aseptic techniques. Insert a catheter introducer of the correct size and then insert the catheter.
- Confirm that the thumbknob is pulled back completely before insertion. Under fluoroscopic guidance, advance the catheter until the tip is in the intended cardiac structure. The catheter tip can be deflected to facilitate positioning by using the proximal control mechanism to vary tip curvature. Pushing the thumbknob forward causes the catheter tip to bend (curve); when the knob is pulled back, the tip straightens.

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Prior to removal of the catheter, confirm that the thumbknob has been pulled back completely. Remove the catheter and dispose of it in an appropriate manner. Do not resterilize and reuse.

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WARNINGS

Cardiac catheterization procedures which generate electrograms and use fluoroscopy in aiding catheter positioning present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging. Cardiac catheterization should only be performed after adequate attention has been given to potential radiation exposure associated with the procedure, and steps taken to minimize this exposure. Careful consideration must therefore be given for the use of this catheter in pregnant women.

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PRECAUTIONS

- Do not attempt to operate the NOGASTAR Catheter prior to completely reading and understanding these instructions for use.
- To avoid thromboemboli, intravenous heparin is used when entering the left side of the heart during cardiac catheterization procedures.
- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped laboratory.
- When using the NOGASTAR Catheter with conventional systems (using fluoroscopy to determine catheter tip location), or with the NOGA System, careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- Always pull the thumbknob of the catheter back before insertion or withdrawal to assure that the catheter tip assumes its original shape.
- Store in a cool, dry, dark place.
- The sterile packaging and catheter should be inspected prior to use.
- The NOGASTAR Catheter is intended for single patient use only.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or crossinfection, including, but not limited to, the transmission of

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infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

- Make sure the small vent hole(s) at the connector end of the handpiece are patent; if clogged, piston could force air into the catheter shaft and possibly into the bloodstream.
- Regularly inspect and test re-usable cables and accessories.
- Federal (USA) law restricts this device to sale by or on the order of a physician.

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There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the product(s) described herein. Under no circumstances shall Biosense Webster,

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Inc. or its affiliated companies, be liable for any special, direct, incidental, consequential or other damages other than as expressly provided by specific law.

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Medical device is stored in a cool, dark, dry place at a room temperature, at a place protected from direct sunlight at temperature

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from $+5^{\circ}\text{C}$ до $+25^{\circ}\text{C}$ and air relative humidity from 20% to 80% (without condensation).

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STERILIZATION

- STERILE. Sterilized using ethylene oxide.
- For single use only.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or crossinfection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

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WARNING: • To reduce the risk of entrapping cardiac structures in the mapping-electrode portion of the catheter, place the LASSO 2515 Variable Circular Mapping Catheter by torquing (or rotating) the shaft in a clockwise motion only.

- When not in regions intended for mapping, manipulate the catheter with the loop in the fully expanded (i.e. 25 mm diameter) position to further reduce the risk of entrapping cardiac structures.
- **STERILE.** Sterilized with ethylene oxide gas.
- For one single use only.
- Do not autoclave.
- Do not use if the package is open or damaged.

DESCRIPTION

The LASSO 2515 Variable Circular Mapping Catheter-Electrode diagnostic (further – electrode) has been designed to facilitate electrophysiological mapping of the atria of the heart. It is deployed in the right or left atrium through an 8F guiding sheath, similar to Biosense Webster's PREFACE Guiding Sheaths. This deflectable catheter consists of a 4F circular spine on its distal tip, with platinum electrodes that can be used for stimulation and recording. The LASSO 2515 Variable Circular Mapping Catheter features a Nitinol loop design that allows the expansion and contraction of the loop to custom-fit veins with different sizes, ranging from 25mm to 15mm diameter ($\pm 15\%$).

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

With the catheter pointing away from you, rotate the handle to the right to reduce the diameter of the loop. Rotating the handle to the left increases the loop diameter (see Figure 1). When the handle is fully rotated to the right, the minimum 15mm diameter is attained.

Likewise, when the handle is rotated fully to the left, the maximum diameter of 25mm is attained. Pushing forward on the catheter thumbknob deflects the tip of the catheter. When the thumbknob is pulled back, the catheter tip straightens.

The catheter interfaces with standard recording equipment via interface cables with the appropriate connectors. Consult the local distributor or manufacturer for the appropriate interface cables.

Figure 1 – LASSO 2515 Variable Circular Mapping Catheter electrodes.

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INDICATIONS FOR USE

The catheter is indicated for multiple electrode electrophysiological mapping of the cardiac structures of the heart, i.e. recording or stimulation only. The LASSO 2515 Variable Circular Mapping Catheter is designed to obtain electrograms in the atrial region of the heart.

CONTRAINDICATIONS

- The LASSO 2515 Variable Circular Mapping Catheter has not been shown to be safe and effective for radio frequency (RF) ablation.
- Use of the catheter may not be appropriate for use in patients with prosthetic valves. A relative contraindication for cardiac catheter procedures is active systemic infection.
- The transseptal approach is contraindicated in patients with left atrial thrombus or myxoma, or interatrial baffle or patch.
- The retrograde approach is contraindicated because of risk of entrapping the LASSO 2515 Variable Circular Mapping Catheter in the left ventricle or valvular apparatus. The LASSO 2515 Variable Circular Mapping Catheter is not recommended for use in the ventricles.

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INSTRUCTIONS FOR USE

- Remove the catheter from its package and place it in a sterile work area.
- Follow standard practice for vessel puncture, guidewire insertion and guiding sheath use and aspiration per its Instructions for Use.
- Connect the interface connectors to the appropriate recording equipment.

NOTE: Read pacing and recording equipment operator manual for proper set up and operation.

- Confirm that the thumbknob is pulled back completely before insertion and that the loop-contraction mechanism is not activated, ensuring minimal tension to the Nitinol loop.

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- Adjust the loop diameter with the handle grip. By rotating the handle grip to the right, the loop contracts; by rotating the grip to the left, the loop expands (see Figure 1).

- Adjust the radius of curvature as necessary by manipulating the thumbknob. Pushing the thumbknob forward causes the catheter tip to bend (curve); when the knob is pulled back, the tip straightens.
- Prior to removal of the catheter, confirm that the thumbknob has been pulled back completely and that the loop is in a fully relaxed position (handle grip rotated fully to the left). Remove the catheter through the guiding sheath and dispose of it in an appropriate manner. Remove the guiding sheath, vessel dilator and guidewire as a unit per its Instructions for Use. Do not resterilize and reuse.
- If there are any questions regarding the use or performance of this product, please consult with the local distributor or the manufacturer.

WARNINGS

- Cardiac catheterization procedures present the potential for significant x-ray exposure, which can result in acute radiation injury, as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging.
- Cardiac catheterization should only be performed after adequate attention has been

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given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure.

- Careful consideration must be given for the use of this catheter in pregnant women.
- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.
- Do not expose the catheter to organic solvents such as alcohol.
- Do not autoclave the catheter.
- To avoid potential damage to anatomical structures, do not attempt to pull the catheter, or withdraw it into the sheath, with the loop in a contracted position. The loop should be fully relaxed (handle grip rotated fully to the left) to minimize tension applied to the Nitinol structure.
- Do not introduce the LASSO 2515 Variable Circular Mapping Catheter tip folded into the guiding sheath.
- The LASSO 2515 Variable Circular Mapping Catheter is recommended for use with the Biosense Webster PREFACE Braided Guiding Sheath.

Note: Do not use the LASSO 2515 Variable Circular Mapping Catheter in conjunction with transseptal sheaths featuring side holes larger than 1.25 mm in diameter.

PRECAUTIONS

- The catheter should be stored in its original packaging and in a cool, dry, dark place until it is used.

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- Read all accessory operating instructions prior to connection of the LASSO 2515 Variable Circular Mapping Catheter.

- Do not attempt to operate the LASSO 2515 Variable Circular Mapping Catheter prior to completely reading and understanding these Instructions for Use.

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- Careful manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance through a guiding sheath. Do not use excessive force to advance or withdraw the catheter through the guiding sheath, when resistance is encountered. In addition, extra care should be taken while inserting, aspirating and manipulating the guiding sheath.

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- The LASSO 2515 Variable Circular Mapping Catheter is intended for single patient use only.

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- To avoid char formation on the LASSO 2515 Variable Circular Mapping Catheter rings, do not apply RF energy when the ablation catheter is in contact with one or more of the LASSO 2515 Variable Circular Mapping Catheter electrodes.

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There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Product(s) described herein. Under no circumstances shall Biosense Webster, Inc.

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or its affiliated companies, be liable for any special, direct, incidental, consequential or other damages other than as expressly provided by specific law.

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DESCRIPTION

The CELSIUS RMT Diagnostic/Ablation Steerable Tip is a multi-electrode device that can be used for electrophysiological mapping and radio frequency (RF) ablations of the heart.

The catheter incorporates a temperature sensor in the tip electrode that permits temperature monitoring during radiofrequency ablation. The temperature sensor is embedded in the tip electrode and electrical connection to the sensor is made via an accessory cable. The tip electrode also serves to deliver radiofrequency current from the radiofrequency lesion generator to the desired ablation site. The catheter has a high torque proximal shaft with a steerable tip section containing an array of platinum electrodes, and a thermocouple temperature sensing capability. Tip deflection is controlled by the interaction of magnets placed in the catheter tip and shaft, and the magnetic field created by the external Stereotaxis Magnetic Navigation System (MNS).

When used in conjunction with standard

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recording equipment, compatible RF generator and a dispersive pad reference electrode, all four electrodes of the CELSIUS RMT diagnostic/ablation Catheter Steerable Tip are capable of intra-cardiac recording and stimulation, while only the tip electrode can deliver RF energy.

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INDICATIONS FOR USE

The CELSIUS RMT diagnostic/ablation Catheter Steerable Tip is indicated for cardiac electrophysiological mapping (stimulation and recording) and, when used in conjunction with a radiofrequency generator, for cardiac ablation.

CONTRAINDICATIONS

Do not use this device:

- in patients with active systemic infection;
- via the transseptal approach in patients with left atrial thrombus or myxoma, or interatrial baffle or patch;
- via the retrograde transaortic approach in patients with aortic valve replacement.

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A number of serious adverse reactions have been documented for catheter ablation procedures including pulmonary embolism, myocardial infarction, stroke, cardiac tamponade, and death. The following complications have been reported in the literature:

- Catheterization/catheter procedure related:
vascular bleeding/local hematomas, thrombosis, AV fistula, pseudoaneurysm, thromboembolism and vasovagal reactions, cardiac perforation, tamponade, thrombi, air embolism, arrhythmias and valvular damage, pneumothorax and hemothorax.
- Radiofrequency related:
chest pain/discomfort, ventricular tachyarrhythmia, TIA, CVA, complete heart block, coronary artery spasm, coronary artery thrombosis, coronary artery dissection, cardiac thromboembolism, pericarditis, cardiac perforation/tamponade, valvular damage and increased phosphokinase level.

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DIRECTIONS FOR USE

Physician Training

Physicians must be familiar with the techniques and appropriately trained for cardiac

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Detailed Device Description

The CELSIUS RMT Diagnostic/Ablation Catheter Steerable Tip is available in a 7 French shaft size with a single 4 mm platinum tip-electrode and three platinum ring electrodes. Each of the ring electrodes are spaced from the tip electrode (or 2, 5 & 2mm from each other) and are capable of performing both intracardiac recording and stimulation. However, the tip electrode is the only electrode that can be used to deliver RF energy to the desired tissue site. Conductor wires from the electrodes pass through the catheter handle and terminate in a Redel type 10-pin connector.

The shaft and tip sections of the CELSIUS RMT Diagnostic/Ablation Catheter Steerable Tip are constructed from a Pebax blended material. All standard versions of the catheter have a combined usable shaft insertion length of 125 ± 5 cm.

Handling and Preparation

Using aseptic technique, remove the catheter from its package and place it in a sterile working area. Inspect the catheter carefully for electrode integrity and overall condition.

Use During the Mapping and Ablation

1. Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
2. Connect the catheter to the recording equipment and/or a compatible RF Generator using the appropriate interface cables.

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3. Confirm that the magnetic fields are removed.

4. Insert catheter and advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning.

5. Use the MNS (Magnetic Navigation System) to orient the distal tip of the catheter, while advancing and retracting the catheter.

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7. RF current may be re-applied to the same or alternate sites using the same catheter.

8. Prior to removal of the catheter, confirm that the magnetic fields are removed.

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The CELSIUS RMT Diagnostic/Ablation Catheter Steerable Tip should be used only with an RF generator which has been shown to be safe and effective for cardiac ablation. Refer to the RF generator manual for detailed generator operating instructions for RF catheter ablation.

If radiofrequency current is interrupted for either a temperature or an impedance rise (a preset limit is exceeded), the catheter should be removed, and the tip cleaned of coagulum. When cleaning the tip electrode, be careful not to twist the tip electrode with respect to the catheter shaft; twisting may damage the tip electrode bond and loosen the tip electrode.

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Use the appropriate Biosense Webster accessory cable to connect the CELSIUS RMT Diagnostic/Ablation Catheter Steerable Tip to a compatible RF generator.

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FIGURE 1. Schematic of the CELSIUS RMT Diagnostic/Ablation Catheters Steerable Tip, Interface Cables and Compatible RF Generators

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FIGURE 2. Schematic of the CELSIUS RMT Diagnostic/Ablation Catheters Steerable Tip, Interface Cables and Compatible RF Generators

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FIGURE 3. Schematic of the CELSIUS RMT Diagnostic/Ablation Catheters Steerable Tip, Interface Cables and Compatible RF Generators

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WARNINGS

- Significant x-ray exposure can result in acute radiation injury as well as dose-related risk for somatic and genetic effects

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(dose = duration of the fluoroscopic imaging X x-ray beam intensity). Take all appropriate measures to minimize x-ray exposure to both patients and clinical staff. Careful consideration should be given to the use of this device in pregnant women.

- **Ablation from within a coronary artery can cause myocardial injury and death.** Adequate fluoroscopic visualization is necessary during the transaortic approach to avoid placement of the catheter in the coronary vasculature.

- **Stroke or myocardial infarction** may occur in patients undergoing left-sided ablation procedures. Patients should be closely monitored during the postablation period for clinical manifestations of embolic events.

- **Implantable pacemakers and implantable cardioverter/defibrillators (ICDs)** may be adversely affected by RF ablation. ICDs should be deactivated during ablation. Have temporary external sources of pacing and defibrillation available during ablation. Exercise extreme caution during ablation when in close proximity to device leads and perform a complete analysis of the implanted device function after ablation.

- Complete AV block can occur when ablating septal accessory pathways or in the treatment of AVNRT. Closely monitor AV conduction during RF energy delivery

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and immediately terminate energy delivery if partial or complete AV block is observed. Using catheters with distal pair electrode spacing greater than 2 mm may increase the risk of AV nodal damage.

- To reduce the risk of brachial plexus injury, physicians are advised to position the patient's arms inferiorly in the normal position with the hands down by the hips as opposed to superiorly.
- Follow all warnings and precautions applicable to the Stereotaxis Magnetic Navigation System, especially those regarding magnetic objects in the catheter laboratory. Refer to Stereotaxis Magnetic Navigation System Operations Manual for user precautions.

PRECAUTIONS

- Cardiac ablation procedures should be performed only by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- Do not attempt to operate the CELSIUS RMT Diagnostic/Ablation Catheter Steerable Tip or the RF generator before completely reading and understanding the applicable directions for use.
- The long-term risks of protracted fluoroscopy have not been established. Careful

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consideration must therefore be given for the use of the device in prepubescent children.

- The long-term risks of lesions created by RF ablation have not been established. In particular, any long-term effects of lesions in proximity to the specialized conduction system or coronary vasculature are unknown.
- When radiofrequency current is interrupted for either a temperature or an impedance rise (the set limit is exceeded), the catheter should be removed, and the tip cleaned of coagulum. When cleaning the tip electrode, be careful not to twist the tip electrode with respect to the catheter shaft; twisting may damage the tip electrode bond and loosen the tip electrode.
- Apparent low power output, high impedance reading or failure of the equipment to function correctly at normal settings may indicate faulty application of the dispersive electrode(s) or failure of an electrical lead. Do not increase power before checking for obvious defects or misapplication.

- Do not attempt to operate the CELSIUS RMT Diagnostic/Ablation Catheter Steerable Tip or radiofrequency generator prior to completely reading and understanding the applicable instructions for use.

Catheter Compatibility

- The CELSIUS RMT Diagnostic/Ablation Catheter Steerable Tip is intended for use with a compatible RF generator and Biosense Webster accessories only.
- Read and follow the dispersive electrode manufacturer's instructions for use; the use of dispersive electrodes, which meet or exceed ANSI/AAMI requirements (HF18), is recommended.

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Handling and Sterilization

- Store in a cool, dry place.
- SINGLE USE ONLY. Observe "Use By" Date. Sterilized with ethylene oxide gas.

- The sterile packaging and catheter should be inspected prior to use. If the package or the catheter appears damaged, do not use. Contact your local Biosense Webster representative.

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- Catheter damage may occur due to
 - autoclaving
 - resterilizing
 - exposure to organic solvents

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Environmental and EMI

- Catheter materials are not compatible with magnetic resonance imaging (MRI).

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- During ablation procedures, this catheter is used in conjunction with an RF generator. Electromagnetic interference (EMI) produced by the radiofrequency (RF) generator during the delivery of RF power may adversely affect the performance of other equipment.

Precautions during Catheter Use

- The patient should not contact grounded metal surfaces. Use only isolated amplifiers, pacing equipment, and ECG equipment or patient injury or death may occur. Leakage current from any connected device to the patient must not exceed 10 microAmps (uA) under any circumstances.
- Do not use excessive force to advance or withdraw the catheter. Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade.
- Always straighten the catheter tip before insertion or withdrawal of the catheter.
- Do not insert or withdraw the catheter without placing the catheter tip in neutral position.
- Use both fluoroscopy and electrograms to monitor the advancement of the catheter to the area of the endocardium under investigation to avoid vascular or cardiac damage.
- Do not immerse the electrical connector or interface cable in fluids; electrical performance could be compromised.

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• Do not deliver RF energy with catheter outside the target site. The RF generator can deliver significant electrical energy and may cause patient or operator injury.

• Avoid use of electrodes and probes of monitoring and stimulating devices, which could provide paths for high frequency current. Reduce the burn hazard by placing the electrodes and probes as far away as possible from the ablation site and the dispersive electrode.

• In the event of a generator cutoff (impedance or temperature), the catheter must be withdrawn and the tip electrode cleaned of coagulum before RF current is re-applied. Use only sterile saline and gauze pad to clean the tip.

• Do not scrub or twist the tip electrode as damage may cause catheter failure or patient injury.

• Discontinue ablation immediately and replace catheter if tip temperature fails to rise during ablation.

• The catheter measures electrode tip temperature, not tissue temperature. If the generator does not display temperature, verify that the appropriate cable is plugged into the generator. If temperature still is not displayed, there may be a malfunction in the temperature sensing system, which must be corrected before applying RF power.

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**DISCLAIMER OF WARRANTY
AND LIMITATION OF LIABILITY
THERE IS NO EXPRESS OR IM-
PLIED WARRANTY, INCLUDING
WITHOUT LIMITATION ANY IM-
PLIED WARRANTY OF MER-
CHANTABILITY OR FITNESS FOR
A PARTICULAR PURPOSE, ON
THE PRODUCT(S) DESCRIBED
HEREIN. UNDER NO CIRCUM-
STANCES SHALL BIOSENSE WEB-
STER, INC., OR ITS AFFILIATED
COMPANIES, BE LIABLE FOR
ANY SPECIAL, DIRECT, INCI-
DENTAL, CONSEQUENTIAL, OR
OTHER DAMAGES OTHER THAN
AS EXPRESSLY PROVIDED BY
SPECIFIC LAW.
WITHOUT LIMITING THE FORE-
GOING, BIOSENSE WEBSTER,
INC. OR ITS AFFILIATED COMPA-
NIES, SHALL NOT BE LIABLE FOR
ANY SPECIAL, DIRECT, INCI-
DENTAL, CONSEQUENTIAL, OR
OTHER DAMAGES, ARISING OUT
OF THE REUSE OF ANY PROD-
UCT(S) LABELED FOR SINGLE
USE OR WHERE REUSE IS PRO-
HIBITED BY APPLICABLE LAW.**

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STERILIZATION

- STERILE. Sterilized with ethylene oxide gas.
- For single use only.

This device is packaged and sterilized for single use only. Do not reuse, reprocess or re-sterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient

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DO NOT RESTERILIZE.

DESCRIPTION

The WEBSTER Coronary Sinus Catheter-Electrode diagnostic (further – catheter) is a steerable, multi-electrode catheter with a deflectable tip designed to facilitate electrophysiological mapping of the heart. The models with EZ steer are both 7 FR catheters with a usable length of 115 cm. The model without EZ Steer is a 6Fr catheter with a usable length of 115cm. The catheter has a hightorque shaft with a bi-directional deflectable tip section containing an array of platinum electrodes that can be used for stimulation and recording.

Standard features of the EZ steer catheters include a braided bi-directional deflectable tip section with an array of ten platinum electrodes that includes a 2 mm tip dome. Additionally, two asymmetric curve types are available providing two 180° opposed, single plane curves. Currently, the available curves for WEBSTER CS include DF and FJ. A Rocker Lever is used to deflect the tip. The high-torque shaft allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site.

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The WEBSTER CS without EZ Steer has standard features including a braided 6FR deflectable tip section with an array of platinum electrodes that includes a 2 mm tip dome. The braided tip is controlled by a proximal hand piece that features a thumb operated sliding piston, and is offered in various curve types. The high-torque shaft allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site.

Both models with Auto ID are equipped with Electronically Erasable Programmable Read Only Memory (EEPROM) which is used to store unique catheter identification information. Carto® EP Navigation Systems equipped with Auto ID Technology can access the stored information and automatically recognize the catheter information. These catheters can interface with

Carto® EP Navigation Systems equipped with Auto ID Technology via interface cables with the appropriate connectors.

CONSULT THE LOCAL
DISTRIBUTOR OR MANUFACTURER
FOR THE APPROPRIATE INTERFACE
CABLES. REFER TO THE
APPROPRIATE Carto® EP
NAVIGATION SYSTEM USER
MANUAL FOR DETAILS ON SYSTEM
SETTINGS RELATED TO THE
CATHETER.

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INDICATIONS FOR USE

The WEBSTER CS Catheter is indicated for electrophysiological mapping of cardiac structures; i.e., stimulation and recording only. The catheter is designed for use in the coronary sinus.

CONTRAINDICATIONS

The WEBSTER CS Catheter has not been shown to be safe and effective for electrical ablation.

- Use of the catheter may not be appropriate for patients with prosthetic valves.
- A relative contraindication for cardiac catheter procedures is active systemic infection.

Do not use this device:

- Via the transseptal approach in patients with left atrial thrombus or myxoma or interatrial baffle or patch;
- Via the retrograde approach in patients with aortic valve replacement.

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The following complications associated with cardiac catheterization have also been

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reported in literature: vascular bleeding, local hematomas, thrombosis, AV fistula, pseudoaneurysm, thromboembolism, vasovagal reactions, cardiac perforation, tamponade, thrombi, air embolism, arrhythmias, valvular damage, pneumothorax and hemothorax.

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INSTRUCTIONS FOR USE

Physician Training

Physicians must be familiar with the techniques and appropriately trained for cardiac mapping procedures. All mapping procedures must be performed in a fully equipped electrophysiology laboratory.

Compatible Accessories:

Use appropriate Biosense Webster accessory cables to connect the WEBSTER CS Catheter to standard recording equipment for the standard model and to the Patient Interface Unit of the appropriate CARTO® EP Navigation System for the model with Auto ID and connect the catheter to the interface cable.

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Before use, inspect the packaging. Do not use if open or damaged. Using aseptic technique, remove the catheter from its package and place it in a sterile working area. Inspect the catheter carefully for electrode integrity and overall condition.

Operating Instructions

For the models with EZ Steer:

1. Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
2. Connect the catheter to the interface cables and standard recording equipment using the appropriate interface cables or the Patient Interface Unit of the appropriate CARTO® EP Navigation System and connect the catheter to the interface cable.
3. Advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning.
4. Use the Rocker Lever to deflect the catheter tip (Figure 1). When the lever is pulled back from neutral, the tip will deflect relative to the direction of rotation. The amount of deflection is relative to the amount of lever rotation. When the lever is pushed forward, the tip will deflect in the opposite direction. To straighten the tip, return the Rocker Lever to neutral position.



Figure 1

The handle has an adjustable friction control that allows the operator to use the Rocker Lever and deflecting tip in a “free” state or adjust the friction to where the Rocker Lever and tip curve are “locked” in place (Figure 2). This knob is located on the opposite side of the Rocker Lever. Out of the package, the knob will be in the “off”

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
position which allows free movement for the lever and deflecting tip. The amount of friction increases as the Friction Control Knob is rotated clockwise until it reaches the fully “on” position.

Clockwise rotation from the “off” position increases the friction within the Deflection Mechanism

Counter-Clockwise rotation from the “on” position decreases the friction within the Deflection Mechanism

Figure 2

5. Determine that the electrodes are in stable contact with the intended mapping site.

For the model without EZ Steer:

1. Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
2. Connect the interface cable to the Patient Interface Unit of the appropriate CARTO® EP Navigation System and connect the catheter to the interface cable.
3. Confirm that the thumb knob is pulled back completely before insertion (catheter in straight position).
4. Advance the catheter into the coronary sinus. Use both fluoroscopy and electrograms to aid in proper positioning.
5. Adjust the radius of curvature as necessary by manipulating the thumb knob. Pushing the thumb knob forward causes the catheter tip to bend (curve); when the knob is pulled back, the tip straightens.
6. Prior to removal of the catheter, confirm that the thumb knob has been pulled back completely (straightened). Remove the catheter and dispose of it in an appropriate manner.

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WARNINGS

- Take all appropriate measures to minimize x-ray exposure to both patients and clinical staff. Significant x-ray exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects due to the x-ray beam intensity and duration of the fluoroscopic imaging. The long-term risks of protracted fluoroscopy have not been established. Careful consideration must therefore be given for the use of the device in prepubescent children and pregnant women.
- Do not attempt to operate the WEBSTER CS Catheter prior to completely reading and understanding the applicable Instructions for Use.
- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.
- For the models with EZ Steer always place the Rocker Lever in the neutral position to straighten the catheter tip before insertion or withdrawal of the catheter. For the model without EZ Steer; Always ensure that the thumb knob is pulled back completely to straighten the catheter tip before insertion or withdrawal of the catheter.

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PRECAUTIONS

Handling and Sterilization

**FOR SINGLE USE ONLY. DO NOT
RESTERILIZE.**

Catheter damage may occur due to:

- autoclaving
- resterilizing
- exposure to organic solvents
- immersing proximal handle or cable connector in fluids
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

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- Sterilized with ethylene oxide gas.
- The sterile packaging and catheter should be inspected prior to use. If the package or the catheter appears damaged, do not use. Contact your local Biosense Webster representative.

HOW SUPPLIED

- The available curves for WEBSTER CS with EZSteer Technology Catheters include DF and FJ.
- The available curves for WEBSTER CS Catheters include D, F and J.
- Appropriate interface cables are supplied separately.

Packaging

The WEBSTER CS Catheter is supplied STERILE (EO). The catheter is secured in a two-piece thermoform tray and placed into a Tyvek/Nylon film pouch, sealed, and placed in a box. Both the pouch and thermoform tray are sterile unless the package is damaged or opened.

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- STERILE. Sterilized with ethylene oxide gas.
- For one single use only.
- Do not autoclave.
- Do not use if the package is open or damaged.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician,

DESCRIPTION

The LASSO Deflectable Circular Mapping Catheter-Electrode Diagnostic (further - catheter) have been designed to facilitate electrophysiological mapping of the atria of the heart. It is deployed in the right or left atrium through an 8F guiding sheath, similar to Biosense Webster's PREFACE Guiding Sheaths. The LASSO deflectable catheter consists of a 3F circular spine on its distal tip, with platinum electrodes that can be used for stimulation and recording. The LASSO AUTO ID model also comes in a 4F circular spine model. Models are available with different circular spine diameters. By pushing forward on the catheter thumbknob, the tip is deflected; when the thumbknob is pulled back, the catheter tip straightens. The Auto ID catheter is equipped with Electronically Erasable Programmable Read Only Memory (EEPROM) which is used to store unique catheter identification information. CARTO EP Navigation Systems equipped with Auto ID Technology can access the stored information and automatically recognize the catheter information.

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The catheter interfaces with CARTO EP Navigation Systems equipped with Auto ID Technology via interface cables with the appropriate connectors.

The standard LASSO catheter interfaces with standard recording equipment via interface cables with the appropriate connectors.

REFER TO THE APPROPRIATE CARTO EP NAVIGATION SYSTEM USER MANUAL FOR DETAILS ON SYSTEM SETTINGS RELATED TO THE CATHETER.

The sterile packaging and catheter should be inspected prior to use. If the package or the catheter appears damaged, do not use. Contact your local Biosense Webster representative

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INDICATIONS FOR USE

The catheter is indicated for multiple electrode electrophysiological mapping the cardiac structures of the heart, i.e. recording or stimulation only. The LASSO Circular Mapping Catheter is designed to obtain electrograms in the atrial region of the heart.

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CONTRAINDICATIONS

- The LASSO Circular Mapping Catheter has not been shown to be safe and effective for radio frequency (RF) ablation.
- Use of the catheter may not be appropriate for use in patients with prosthetic valves. A relative contraindication for cardiac catheter procedures is active systemic infection.
- The transseptal approach is contraindicated in patients with left atrial thrombus or myxoma, or interatrial baffle or patch.
- The retrograde approach is contraindicated because of risk of entrapping the LASSO Catheter in the left ventricle or valvular apparatus. The LASSO Catheter is not recommended for use in the ventricles.

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INSTRUCTIONS FOR USE

- Remove the catheter from its package and place it in a sterile work area.
- Follow standard practice for vessel puncture, guidewire insertion and guiding sheath use and aspiration per its Instructions for Use.
- Connect the interface cable to the appropriate recording equipment for the standard model and the Patient Interface Unit of the appropriate CARTO EP Navigation System for the AUTO ID model and connect the catheter to the interface cable.
- Confirm that the thumbknob is pulled back completely before insertion. Advance the catheter through the guiding sheath to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning. Adjust the radius of curvature as necessary by manipulating the thumbknob. Pushing the thumbknob forward causes the catheter tip to bend (curve); when the knob is pulled back, the tip straightens.
- To place LASSO Catheter, (or rotate) shaft in a clockwise motion only.
- Prior to removal of the catheter, confirm that the thumbknob has been pulled back completely. Remove the catheter through the guiding sheath and dispose of it in an appropriate manner. Remove the guiding sheath, vessel dilator and guidewire as a unit per its

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Instructions for Use. Do not resterilize and reuse.

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WARNINGS

- Cardiac catheterization procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging.
- Cardiac catheterization should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure.
- Careful consideration must be given for the use of this catheter in pregnant women.
- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.
- Do not expose the catheter to organic solvents such as alcohol.
- Do not autoclave the catheter.
- Do not introduce the LASSO Catheter tip folded into the guiding sheath.
- The LASSO Deflectable Circular Mapping Catheter is recommended for use with

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the Biosense Webster PREFACE Braided Guiding Sheath.

Note: Do not use the LASSO Deflectable Circular Mapping Catheter in conjunction with transseptal sheaths featuring side holes larger than 1.25mm in diameter.

PRECAUTIONS

- The catheter should be stored in its original packaging and in a cool, dry place until it is used.
- Do not attempt to operate the LASSO Circular Mapping Catheter prior to completely reading and understanding these Instructions for Use.
- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- Careful manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance through a guiding sheath. Do not use excessive force to advance or withdraw the catheter through the guiding sheath when resistance is encountered. In addition, extra care should be taken while inserting, aspirating and manipulating the guiding sheath.
- The sterile packaging and catheter should be inspected prior to use.
- Do not resterilize and reuse. It is intended for single patient use only.

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- To place LASSO Catheter, torque (or rotate) shaft in a clockwise motion only.

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STERILIZATION

- STERILE. Sterilized with ethylene oxide gas.
- For one single use only.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

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SINGLE USE ONLY. DO NOT RESTERILIZE.

DESCRIPTION

The NAVISTAR RMT Diagnostic/Ablation Steerable Tip Catheter-Electrode (further – catheter) is a steerable, multi-electrode catheter designed to facilitate electrophysiological mapping of the heart and to transmit radiofrequency current to the catheter tip electrode for ablation purposes. For mapping, the catheter is used with a RMT compatible CARTO EP Navigation System (a magnetic field location technology) and an appropriate location reference device. For ablation, the catheter is used in conjunction with a compatible radiofrequency (RF) generator and a grounding (dispersive) pad.

The device has a high-torque nylon shaft with a steerable tip section containing an array of platinum electrodes, and a thermocouple temperature sensing capability. Tip deflection is controlled by the interaction of magnets placed in the catheter tip and shaft, and the magnetic field created by the external Stereotaxis Magnetic Navigation System (MNS). All electrodes may be used for recording and stimulation, but only the tip electrode may be used to deliver RF energy from the generator. A lo-

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The catheter interfaces with a RMT compatible CARTO EP Navigation System and a compatible RF generator via interface cables and a Patient Interface Unit (PIU).

For use in mapping procedures, refer to the instructions for the RMT compatible CARTO EP Navigation System.

For use in RF ablation procedures, refer to the RF generator operating instructions.

For use in mapping and ablation procedures, refer to the instructions for the Stereotaxis Magnetic Navigation System.

INDICATIONS AND USAGE

The NAVISTAR RMT Diagnostic/Ablation Steerable Tip Catheter, and related accessory devices are indicated for catheter-based atrial and ventricular electrophysiological mapping, and for use with a compatible RF generator in adults and children four (4) years of age and older in

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creating endocardial lesions during cardiac ablation procedures to treat arrhythmias.

When used with a RMT compatible CARTO EP Navigation System, the NAVISTAR RMT Diagnostic/Ablation Tip Catheter provides location information.

CONTRAINDICATIONS

Do not use this device:

- in patients with active systemic infection;
- via the transseptal approach in patients with left atrial thrombus or myxoma, or interatrial baffle or patch;
- via the retrograde approach in patients with aortic valve replacement.

ADVERSE REACTIONS

A number of serious adverse reactions have been documented for catheter ablation procedures including pulmonary embolism, myocardial infarction, stroke, cardiac tamponade, and death.

The following complications were also noted to have occurred during prior studies or have been reported in the literature:

- Catheterization/catheter procedure related: vascular bleeding/local hematomas, thrombosis, arteriovenous (AV) fistula, pseudoaneurysm, thromboembolism and vasovagal reactions, cardiac perforation, pericardial effusion/tamponade, thrombi, air embolism, arrhythmias and valvular

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damage, pneumothorax and hemothorax, pulmonary edema, hypoxia, pleural effusion, acute respiratory distress syndrome (ARDS), congestive heart failure, aspiration pneumonia, pneumonia, asthmatic attack, hypotension, implantable cardioverter-defibrillator (ICD) lead malfunction, anemia, thrombocytopenia, disseminated intravascular coagulation, epistaxis, systemic infection, urinary tract infection, sedation induced apnea, sedation induced CO₂ retention with lethargy and cholecystitis.

- RF related: chest pain/discomfort, ventricular tachyarrhythmia, transient ischemic attack (TIA), cerebrovascular accident (CVA), complete heart block, coronary artery spasm, coronary artery thrombosis, coronary artery dissection, cardiac thromboembolism, pericarditis, cardiac perforation/tamponade, valvular damage and increased phosphokinase level.
- Unrelated to device or procedure: urinary retention, transient extremity numbness, Parkinson's disease and gastrointestinal diverticulosis.

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INSTRUCTIONS FOR USE

Physician Training

Physicians must be familiar with the techniques and appropriately trained for cardiac mapping and ablation procedures. All mapping and ablation procedures must be performed in a fully equipped electrophysiology laboratory.

Compatible RF Generators and Accessories

The NAVISTAR RMT Diagnostic/Ablation Steerable Tip Catheter should be used only with a legally marketed RF generator, which has been shown to be safe and effective for cardiac ablation. Refer to the RF generator operator's manual for detailed operating instructions for RF ablation including information in choosing between temperature or power control modes.

Specifications for Compatible RF Generators: Refer to the RF generator manual for detailed generator operating instructions for RF catheter ablation.

Accessories: Use appropriate Biosense Webster accessory cables to connect the NAVISTAR RMT Diagnostic/Ablation Steerable Tip Catheter to the RMT compatible CARTO EP Navigation System and a compatible RF generator. Use an appropriate location reference device when using the catheter with the RMT compatible CARTO EP Navigation System. Refer to the CARTO EP Navigation System User's Manual for detailed instructions for information about catheter interface cable, location reference device and mapping procedures.

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Handling and Preparation

Before use, inspect the packaging. Do not use if open or damaged.

Using aseptic technique, remove the catheter from its package and place it in a sterile working area. Inspect the catheter carefully for electrode integrity and overall condition. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.

Operating Instructions

1. Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
2. Connect the catheter to the interface cables, the RMT compatible CARTO EP Navigation System, a compatible RF generator, and standard recording equipment using the appropriate interface cables.
3. Confirm that the magnetic fields are removed.
4. Insert catheter and advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning.
5. Use the MNS (Magnetic Navigation System) to orient the distal tip of the catheter, while advancing and retracting the catheter.

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6. The catheter tip can be deflected to facilitate positioning by using the MNS to vary tip curvature. When it has been determined that the tip electrode is in stable contact with the intended ablation site, the catheter tip electrode connection must be switched from the recording equipment to the RF generator in preparation for delivery of RF current.

7. Determine that the tip electrode is in stable contact with the intended ablation site. Circuit impedance should be approximately 100 ohms upon initiation of RF current.

8. RF current may be re-applied to the same or alternate sites using the same catheter.

9. Prior to removal of the catheter, confirm that the magnetic fields are removed.

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WARNINGS/ PRECAUTIONS

- **Significant x-ray exposure can result in acute radiation injury** as well as dose-related risk for somatic and genetic effects (dose = duration of the fluoroscopic imaging X x-ray beam intensity). Take all appropriate measures to minimize x-ray exposure to both patients and clinical staff. Careful consideration should be given to the use of this device in pregnant women.

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- **Ablation from within a coronary artery can cause myocardial injury and death.** Adequate fluoroscopic visualization is necessary during the transaortic approach to avoid placement of the ablation catheter in the coronary vasculature.

- **Stroke or myocardial infarction** may occur in patients undergoing **left-sided ablation procedures.** Patients should be closely monitored during the postablation period for clinical manifestations of embolic events.

- **Implantable pacemakers and implantable cardioverter/defibrillators (ICDs)** may be adversely affected by RF ablation. ICDs should be deactivated during ablation. Have temporary external sources of pacing and defibrillation available during ablation. Exercise extreme caution during ablation when in close proximity to device leads, and perform a complete analysis of the implanted device function after ablation.

- **Complete AV block** can occur when ablating septal accessory pathways or in the treatment of AVNRT. Closely monitor AV conduction during RF energy delivery and immediately terminate energy delivery if partial or complete AV block is observed. Using catheters with **distal pair electrode spacing greater than 2 mm** may increase the risk of AV nodal damage.

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- **To reduce the risk of brachial plexus injury**, physicians are advised to position the patient's arms inferiorly in the normal position with the hands down by the hips as opposed to superiorly.

- Follow all warnings and precautions applicable to the Stereotaxis Magnetic Navigation System, CARTO System and RF generator especially those regarding magnetic objects in the catheter laboratory. Refer to Stereotaxis Magnetic Navigation System User's Manual for user precautions.

PRECAUTIONS

- Cardiac ablation procedures should be performed only by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- Do not attempt to operate the NAVISTAR RMT Diagnostic/Ablation Steerable Tip Catheter or the RF generator prior to completely reading and understanding the applicable Instructions for Use.
- The long-term risks of protracted fluoroscopy have not been established. Careful consideration must therefore be given for the use of the device in prepubescent children.

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- The long-term risks of lesions created by RF ablation have not been established. In particular, any long-term effects of lesions in proximity to the specialized conduction system or coronary vasculature are unknown.

- When radiofrequency current is interrupted for either a temperature or an impedance rise (the set limit is exceeded), the catheter should be removed, and the tip cleaned of coagulum. When cleaning the tip electrode, be careful not to twist the tip electrode with respect to the catheter shaft; twisting may damage the tip electrode bond and loosen the tip electrode.

- Apparent low power output, high impedance reading or failure of the equipment to function correctly at normal settings may indicate faulty application of the dispersive electrode(s) or failure of an electrical lead. Do not increase power before checking for obvious defects or misapplication.

Catheter Compatibility

- The NAVISTAR RMT Diagnostic/Ablation Steerable Tip Catheter is intended for use with a compatible RF generator and Biosense Webster accessories only in the Niobe magnetic environment.

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SINGLE USE ONLY. DO NOT RE-STERILIZE.

Observe "Use By" Date. Sterilized with ethylene oxide gas.

Catheter damage may occur due to:

- Autoclaving
- Resterilizing
- exposure to organic solvents
- immersing proximal handle or cable connector in fluids

The sterile packaging and catheter should be inspected prior to use. If the package or the catheter appears damaged, do not use. Contact your local Biosense Webster representative.

This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

Environmental and EMI Precautions

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- During ablation procedures, this catheter is used in conjunction with a RF generator. Electromagnetic interference (EMI) produced by the RF generator during the delivery of RF power may adversely affect the performance of other equipment. This device is a non-ionizing emitter of radiation and may cause electromagnetic interference with other devices. In order to minimize the electromagnetic interference, the generator should be placed at least 1 meter away from any other electronic device.

Precautions During Catheter Use

- The patient should not contact grounded metal surfaces. Leakage current from any connected device to the patient must not exceed 10 microAmps (μA) under any circumstances. Use only isolated amplifiers, pacing equipment, and ECG equipment or patient injury or death may occur.

- Do not use excessive force to advance or withdraw the catheter. Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade.

- Always straighten the catheter tip before insertion or withdrawal of the catheter.

- Do not insert or withdraw the catheter without placing the catheter tip in neutral position.

- Use both fluoroscopy and electrograms to monitor the advancement of the catheter to the area of the endocardium under investigation to avoid vascular or cardiac damage.

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- Do not immerse the electrical connector or interface cable in fluids; electrical performance could be compromised.

Precautions During Ablation

- Do not increase power before checking for lead connection and appropriate dispersive electrode application. Effective contact between the patient and the dispersive electrode must be verified whenever the patient is repositioned.

- Do not deliver RF energy with catheter outside the target site. The RF generator can deliver significant electrical energy and may cause patient or operator injury.

- Avoid use of electrodes and probes of monitoring and stimulating devices, which could provide paths for high frequency current. Reduce the burn hazard by placing the electrodes and probes as far away as possible from the ablation site and the dispersive electrode.

- In the event of a generator cutoff (impedance or temperature), the catheter must be withdrawn and the tip electrode cleaned of coagulum before RF current is reapplied. Use only sterile saline and gauze pad to clean the tip.

- Do not scrub or twist the tip electrode as damage may cause catheter failure or patient injury.

- Discontinue ablation immediately and replace catheter if tip temperature fails to rise during ablation.

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- The temperature sensor of the catheter measures electrode tip temperature, not tissue temperature. If the generator does not display temperature, verify that the appropriate cable is plugged into the generator. If temperature still is not displayed, there may be a malfunction in the temperature sensing system, which must be corrected prior to applying RF power.

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PATIENT TREATMENT

Antiplatelet or Anticoagulation Use

To avoid thromboemboli, intravenous heparin is used when entering the left heart during ablation, and many physicians prescribe aspirin, less often warfarin, for about 3 months afterward. No consensus yet exists about the need for short-term anticoagulation after ablation.

Specific Patient Populations

The safety and effectiveness of cardiac ablation has not been established in:

- Asymptomatic patients;
- patients who are pregnant; or
- nursing mothers.

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
PATIENT COUNSELING INFORMATION

Patients may require anticoagulation and/or antiplatelet therapy for an indefinite period based on the patient's condition.

HOW SUPPLIED

The NAVISTAR RMT Diagnostic/Ablation Steerable Tip Catheter is provided sterile (EtO). The device is a 7 F catheter with a usable length of 125 cm, and the following features:

Tip electrode: 4 mm
Connector type: Hypertronics interlocking
Temperature sensor: Thermocouple
French size: 7

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Packaging

The catheter is packaged in a heat-sealed tray. A foam insert is used to secure the catheter handle in place inside the tray package. The sealed tray is then placed into a pouch, forming a double sterile barrier. The sealed pouch

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is placed into an outer folding carton with an instruction for use (IFU) insert.

Storage

The NAVISTAR RMT Diagnostic/Ablation Steerable Tip Catheter must be stored in a cool, dry place. Storage temperature should be from 5 to 25 C (from 41 to 77F).

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This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may

create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

SHELF LIFE

Product and package testing have been conducted to support the "Use By" date printed on the product labels. **DO NOT USE** after the "Use By" date.

Sterility maintaining term - 3 years.

DISPOSAL

Recycle components, or dispose of the product and its residual elements or waste items in accordance with local laws and regulations.

The products have come into contact with

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blood and / or body fluids, disposal and / or destruction according to SanPiN 2.1.7.2790-10 (B class).

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Instruction for Use Catheter-Electrode Navigation Variable Catheter

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

WARNING: To reduce the risk of entrapping cardiac structures in the mapping electrode portion of the catheter, place the LASSO 2515 NAV Variable Catheter by torquing (or rotating) the shaft in a clockwise motion only.

WARNING: When not in regions intended for mapping, manipulate the catheter with the loop in the fully expanded (i.e. 25 mm diameter) position to further reduce the risk of entrapping cardiac structures.

- **STERILE.** Sterilized with ethylene oxide gas.
- For single use only.
- Do not autoclave.
- Do not use if the package is open or damaged.

DESCRIPTION

The LASSO 2515 NAV Variable Catheter-Electrode Navigation (further – catheter) has been designed to facilitate electrophysiological mapping of the atria of the heart with the CARTO 3 EP Navigation System and a reference device. It is deployed in the right or left atrium through an 8F guiding sheath. This deflectable catheter consists of a 4F circular spine on its distal tip, with platinum electrodes that can be used for stimulation and recording.

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The LASSO 2515 NAV Variable Catheter features a Nitinol loop design that allows the expansion and contraction of the loop to custom-fit veins with different sizes, ranging from 25mm to 15mm diameter ($\pm 15\%$).

With the catheter pointing away from you, rotate the handle clockwise to reduce the diameter of the loop. Rotating the handle counter-clockwise increases the loop diameter (see Figure 1). When the handle is fully rotated clockwise, the minimum 15mm diameter is attained. Likewise, when the handle is fully rotated counter-clockwise, the maximum diameter of 25mm is attained. Pushing forward on the catheter thumbknob deflects the tip of the catheter. When the thumbknob is pulled back, the catheter tip straightens.

Figure 1: LASSO 2515 NAV Variable Catheter

The catheter interfaces with standard recording equipment and CARTO 3 EP Navigation System equipment via interface cables with the appropriate connectors. Consult the local distributor or manufacturer for the appropriate interface cables.

For use in mapping procedures, refer to the instructions for the CARTO 3 EP Navigation System.

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For further description of the operation of the CARTO 3 EP Navigation System, refer to the operating instructions for this instrument.

INDICATIONS FOR USE

The LASSO 2515 NAV Variable Catheter is indicated for multiple electrode electrophysiological mapping of the cardiac structures of the heart, i.e. recording or stimulation only. The Catheter is designed to obtain electrograms in the atrial regions of the heart.

The LASSO 2515 NAV Variable Catheter provides location information when used with compatible CARTO EP Navigation Systems version 1.0 or higher.

CONTRAINDICATIONS

- The LASSO 2515 NAV Variable Catheter has not been shown to be safe and effective for radio frequency (RF) ablation.
- Use of the catheter may not be appropriate for use in patients with prosthetic valves. A relative contraindication for cardiac catheter procedures is active systemic infection.
- The transseptal approach is contraindicated in patients with left atrial thrombus or myxoma, or interatrial baffle or patch.
- The retrograde approach is contraindicated because of risk of entrapping the LASSO 2515 NAV Variable Catheter in the left ventricle or valvular apparatus. The LASSO 2515 NAV Variable Catheter is not recommended for use in the ventricles.

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ADVERSE REACTIONS

A number of serious adverse reactions have been documented for cardiac catheterization procedures including pulmonary embolism, myocardial infarction, stroke, cardiac tamponade, and death.

The following complications associated with cardiac catheterization have also been reported in the literature: vascular bleeding, local hematomas, thrombosis, AV fistula, pseudoaneurysm, thromboembolism, vasovagal reactions, cardiac perforation, air embolism, arrhythmias, valvular damage, pneumothorax and hemothorax.

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INSTRUCTIONS FOR USE

- Remove the catheter from its package and place it in a sterile work area.
- Follow standard practice for vessel puncture, guidewire insertion and guiding sheath use and aspiration per its Instructions for Use.

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WARNINGS

1. Cardiac catheterization procedures present the potential for significant x-ray exposure, which can result in acute radiation injury, as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging.
2. Cardiac catheterization should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure.
3. Careful consideration must be given for the use of this catheter in pregnant women.
4. Do not immerse the proximal handle, pig tail or cable connector in fluids; electrical performance could be affected.
5. Do not expose the catheter to organic solvents such as alcohol.
6. Do not autoclave the catheter.
7. To avoid potential damage to anatomical structures, do not attempt to pull the catheter, or withdraw it into the sheath, with the loop in a contracted position. The loop should be fully relaxed (handle grip rotated fully to the left) to minimize tension applied to the Nitinol structure.
8. Do not introduce the LASSO 2515 NAV Variable Catheter tip folded into the guiding sheath.
9. The LASSO 2515 NAV Variable Catheter is recommended for use with an 8F guiding sheath. Note: Do not use the LASSO 2515 NAV Variable Catheter in conjunction with transseptal sheaths featuring side holes larger than 1.25mm in diameter.

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PRECAUTIONS

1. The standard transseptal procedure should be followed during mapping when moving the catheter from the right atrium to the left atrium.
2. To prevent entanglement of the catheter with the valves and to prevent slippage of the catheter into the ventricles, care should be taken when using the catheter in or around the atrioventricular valve region.
3. Prior to mapping, the physician should ensure that the intracardiac signals are recorded by all of the electrodes on the loop of the catheter.
4. Do not use in the proximity of magnetic resonance imaging (MRI) equipment because the MRI equipment may induce movement of the catheter, resulting in perforation.
5. The catheter should be stored in its original packaging and in a cool, dry place until it is used.
6. Read all accessory operating instructions prior to connection of the LASSO 2515 NAV Variable Catheter.
7. Do not attempt to operate the LASSO 2515 NAV Variable Catheter prior to completely reading and understanding these Instructions for Use.
8. Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.
9. Careful manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance through a guiding sheath. Do not use excessive force to advance or withdraw the catheter through the guiding

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sheath, when resistance is encountered. In addition, extra care should be taken while inserting, aspirating and manipulating the guiding sheath.

10. The sterile packaging and catheter should be inspected prior to use.

11. The LASSO 2515 NAV Variable Catheter is intended for single patient use only.

12. This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

13. Do not resterilize and reuse.

14. Always pull the thumbknob of the catheter back before insertion or withdrawal to assure that the catheter tip assumes its original shape.

15. Do not apply RF energy when the ablation catheter is in contact with one or more of the LASSO 2515 NAV Variable Catheter electrodes.

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Descriptions or specifications appearing in Biosense Webster, Inc. printed matter, including this publication, are informational only and meant solely to generally describe the product at the time of manufacture and are not made or given

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as a warranty of the prescribed product in any way.

CONDITIONS OF TRANSPORTATION, STORAGE AND OPERATION

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Medical device is stored at a room temperature, at a place protected from direct sunlight at temperature from +5 to +25 °C and air relative humidity from 20% to 80% (without condensation).

The medical device is transported in cool, dry, dark place, at temperature from -40 °C to +54 °C and air relative humidity from 20% to 80% (without condensation).

Medical device in operation are stable against influence of climatic factors for temperature condition from +32 °C to +42 °C and to influence of biological liquids.

STERILIZATION

- STERILE. Sterilized with ethylene oxide gas.
- For single use only.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or

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Instruction for Use Catheter-Electrode Navigation

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WARNING: To reduce the risk of entrapping cardiac structures in the mapping electrode portion of the catheter, place the LASSO NAV Catheter by torquing (or rotating) the shaft in a clockwise motion only.

- STERILE. Sterilized with ethylene oxide gas.
- For single use only.
- Do not autoclave.
- Do not use if the package is open or damaged.

DESCRIPTION

The LASSO NAV Catheter-Electrode Navigation (further – catheter) has been designed to facilitate electrophysiological mapping of the atria of the heart. It is deployed in the right or left atrium through an 8 Fr guiding sheath. This deflectable catheter consists of a 4.5 Fr circular spine on its distal tip, with platinum electrodes that can be used for stimulation and recording. Models are available with different circular spine diameters.

The LASSO NAV Catheter provides location information when used with a compatible CARTO EP Navigation System.

The catheter interfaces with standard recording equipment and CARTO EP Navigation System equipment via interface ca-

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bles with the appropriate connectors. Consult local distributor or manufacturer for the

appropriate interface cables. For use in mapping with a CARTO EP Navigation System, refer to the instructions for the CARTO EP Navigation System.

For further description of the operation of the CARTO EP Navigation System, refer to the operating instructions for this instrument.

INDICATIONS FOR USE

The LASSO NAV Catheter is indicated for multiple electrode electrophysiological mapping of the cardiac structures of the heart, i.e. recording or stimulation only. This Catheter is designed to obtain electrograms in the atrial regions of the heart.

The LASSO NAV Catheter provides location information when used with a compatible CARTO EP Navigation System.

CONTRAINDICATIONS

- The LASSO NAV Catheter has not been shown to be safe and effective for radiofrequency (RF) ablation.
- Use of the catheter may not be appropriate for use in patients with prosthetic valves.
- A relative contraindication for cardiac catheter procedures is active systemic infection.

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- The transseptal approach is contraindicated in patients with left atrial thrombus or myxoma, or interatrial baffle or patch.

- The retrograde approach is contraindicated because of risk of entrapping the LASSO NAV Catheter in the left ventricle or valvular apparatus.

- The LASSO NAV Catheter is contraindicated for use in the ventricles.

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USAGE

- Remove the catheter from its package and place it in a sterile work area.

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- The catheter is recommended for use with an 8 Fr atraumatic soft tipped sheath. Confirm compatibility of the catheter with the guiding sheath by fully inserting and withdrawing the catheter through the irrigated sheath before clinical use. If excessive force is required or interference between the catheter and sheath is observed, use an alternate guiding sheath to avoid damaging the catheter or sheath.

- Follow standard practice for vessel puncture, guidewire insertion and guiding sheath use and aspiration per its Instructions for Use.

- Connect the interface connectors to the appropriate recording equipment.

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- Confirm that the thumbknob is pulled back completely before insertion ensuring that the tip is not deflected (straight position).

- To place LASSO NAV Catheter, torque (or rotate) shaft in a clockwise motion only.

- Pushing the thumbknob forward causes the catheter tip to bend (curve); when the knob is pulled back, the tip straightens.

- Prior to removal of the catheter, confirm that the thumbknob has been pulled back completely (tip straightened). Remove the catheter through the guiding sheath and dispose of it in an appropriate manner. Remove the guiding sheath, vessel dilator and guidewire as a unit per its Instructions for Use. Do not resterilize and

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WARNINGS

- Cardiac catheterization procedures present the potential for significant x-ray exposure, which can result in acute radiation injury, as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging.
- Cardiac catheterization should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure.
- Careful consideration must be given for the use of this catheter in pregnant women.
- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.
- Do not expose the catheter to organic solvents such as alcohol.
- Do not autoclave the catheter.
- Do not introduce the LASSO NAV Catheter tip folded into the guiding

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sheath.

- The LASSO NAV Catheter is recommended for use with an 8 Fr atraumatic soft tipped guiding sheath.

Note: Do not use the LASSO NAV Catheter in conjunction with transseptal sheaths featuring side holes larger than 1.25 mm in diameter.

PRECAUTIONS

- The standard transseptal procedure should be followed during mapping when moving the catheter from the right atrium to the left atrium.

- To prevent entanglement of the catheter with the valves and to prevent slippage of the catheter into the ventricles, care should be taken when using the catheter in or around the atrial-ventricle valve region.

- Prior to use the physician should ensure that the intracardiac signals are recorded by all of the electrodes on the loop of the catheter.

- Do not use in the proximity of magnetic resonance imaging (MRI) equipment because the MRI equipment may induce movement of the catheter that could result in perforation.

- The catheter should be stored in its original packaging and in a cool, dry place until it is used.

- Read all accessory operating instructions prior to connection of the LASSO NAV Catheter.

- Do not attempt to operate the LASSO NAV Catheter prior to completely reading and understanding these Instructions for Use.

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- To place LASSO NAV Catheter, torque (or rotate) shaft in a clockwise motion only.
- Do not apply RF energy when the ablation catheter is in contact with one or more of the LASSO NAV Catheter electrodes.

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There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on products described here(in). Under no circumstances shall Biosense Webster, Inc. or its affiliated companies, be liable for any special, direct, incidental, consequential or other damages other than as expressly provided by specific law.

Without limiting the foregoing, Biosense Webster, Inc. or its affiliated companies, shall not be liable for any special, direct,

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STERILIZATION

- STERILE. Sterilized with ethylene oxide gas.
- For single use only.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or crossinfection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

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Sterility maintaining term - 3 years.

DISPOSAL

Recycle components, or dispose of the product and its residual elements or waste

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items in accordance with local laws and regulations.

The products have come into contact with blood and / or body fluids, disposal and / or destruction according to SanPiN 2.1.7.2790-10 (B class).

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- STERILE. Sterilized with ethylene oxide (EO) gas.
- For single use only.
- Do not autoclave.
- Do not use if the package is open or damaged.
- CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

DESCRIPTION

The PENTARAY NAV High-Density Mapping Catheter-Electrode navigation is designed to facilitate electrophysiological mapping of the heart with the CARTO 3 EP Navigation System. It is designed for deployment in a heart chamber through an 8F guiding sheath. This deflectable catheter consists of multiple 3F spines on its distal tip, each spine having multiple platinum electrodes that are used for stimulation and recording. A magnetic location sensor embedded in the deflectable tip transmits location information to the CARTO 3 EP Navigation System. The catheter has two electrodes on the deflectable tip to provide for visualization of the tip when used with the CARTO 3 EP Navigation system. Pushing forward on the catheter thumbknob deflects the tip; pulling back on the thumbknob straightens the tip. This device includes an irrigation lumen for connection to a source of continuous drip anticoagulant fluid.

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This catheter interfaces with standard recording equipment and the CARTO 3 EP Navigation System via interface cables with the appropriate connectors.

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INDICATIONS FOR USE

The PENTARAY NAV High-Density Mapping Catheter is indicated for multiple electrode electrophysiological mapping of cardiac structures in the heart, i.e., recording or stimulation only. This catheter is intended to obtain electrograms in the atrial and ventricular regions of the heart. The PENTARAY NAV High-Density Mapping Catheter provides location information when used with compatible CARTO 3 EP Navigation Systems. (This catheter is not compatible with CARTO 3 EP Navigation Systems prior to Version 3.0).

CONTRAINDICATIONS

- The PENTARAY NAV High-Density Mapping Catheter has not been shown to be safe and effective for radio frequency

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(RF) ablation.

- Use of this catheter may not be appropriate for use in patients with prosthetic valves. A relative contraindication for cardiac catheter procedures is active systemic infection.

- The transseptal approach is contraindicated in patients with intracardiac thrombus or myxoma, or interatrial baffle or patch.

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The following complications associated with cardiac catheterization have also been reported in the literature: vascular bleeding, local hematomas, thrombosis, AV fistula, pseudoaneurysm, thromboembolism, vasovagal reactions, cardiac perforation, air embolism, arrhythmias, valvular damage, pneumothorax, and hemothorax.

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USEGE

- Remove the PENTARAY NAV High-Density Mapping Catheter from its

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
package and place it in a sterile work area.


- Follow standard practice for vessel puncture, guidewire insertion, and guiding sheath use and aspiration per its Instructions for Use.
- Connect the interface connectors to the CARTO 3 EP Navigation System.
- Flush the catheter with heparinized saline prior to insertion into the body.
- Always follow standard practices of using a continuous drip of anticoagulant fluid under pressure through the proximal Luer connector when the device is in the body.
- Confirm that the thumbknob is pulled back completely before insertion.
- Advance the insertion tube along the catheter shaft to collapse the spines together prior to insertion into the sheath (see illustration). Do not bend the spines backward. After insertion, slide the insertion tube back toward the handle.

- Advance the catheter through the guiding sheath to the area of the endocardium under evaluation. The spines will naturally return to their expanded pattern as they emerge from the tip of the guiding sheath. Use both fluoroscopy and electrograms to aid in proper positioning. Adjust the radius of curvature as necessary by manipulating the thumbknob. Pushing the thumbknob forward causes the catheter tip

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to curve; pulling the thumbknob back straightens the catheter tip.

- Under fluoroscopy, catheter spine 'A' can be identified by the marker band near ring 2. Catheter spine 'B' can be identified by the marker band near ring 8. Once spines A and B have been identified, all other spines and electrodes may be identified as follows:

- For questions regarding the use or performance of this product, please consult with the local distributor or the manufacturer.

- Prior to removal of the catheter, confirm that the thumbknob has been pulled back completely to straighten the tip. Remove the catheter from the guiding sheath and dispose of it in an appropriate manner.

- Do not resterilize and reuse this device.

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Cardiac catheterization procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging.

Cardiac catheterization should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure.

Careful consideration must be given for the use of this catheter in pregnant women.

Do not immerse the proximal handle, pig tail or cable connector in fluids; electrical performance could be affected.

Do not expose this catheter to organic solvents such as alcohol.

Do not autoclave this catheter.

Flush the catheter with heparinized saline prior to insertion into the body.

Do not introduce the PENTARAY NAV Catheter into a guiding sheath with its distal spines folded backwards (i.e., toward the handle). Collapse the spines together using the insertion tube prior to insertion.

The PENTARAY NAV Catheter is recommended for use with an 8F guiding sheath.

Note: Do not use the PENTARAY NAV Catheter in conjunction with transseptal sheaths featuring side holes larger than 1.25 mm in diameter.

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or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

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- Flush the catheter with heparinized saline prior to insertion into the body. Always follow standard practices of using a continuous drip of anticoagulant fluid under pressure through the proximal Luer connector when the device is in the body.

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There is no express or implied warranty, including without limitation any implied

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warranty of merchantability or fitness for a particular purpose, on the product(s) described herein. Under no circumstances shall Biosense Webster, Inc. or its affiliated companies, be liable for any special, direct, incidental, consequential or other damages other than as expressly provided by specific law.

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STERILIZATION

- STERILE. Sterilized with ethylene oxide (EO) gas.
- For single use only.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

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Catheters-Electrodes Therapeutic

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DESCRIPTION

The Diagnostic/Ablation Catheter-Electrode therapeutic (further – Catheter) is a multi-electrode catheter with a deflectable tip designed to facilitate electrophysiological mapping of the heart and to transmit radiofrequency current to the catheter tip electrode for ablation purposes. For ablation, the catheter is used in conjunction with a radiofrequency generator and a dispersive pad (reference electrode). Two models of the catheter are available: the standard model and the temperature sensing model. The temperature sensing model incorporates a temperature sensor in the tip electrode which permits temperature monitoring during radiofrequency ablation. The temperature sensor is embedded in the large tip electrode and electrical connection to the sensor is made via an accessory cable. The catheter has a high-torque shaft with a deflectable tip section containing an array of platinum electrodes. All of the electrodes may be used for recording and stimulation purposes. The large tip electrode also serves to deliver radiofrequency (RF) current from the radiofrequency lesion generator to the desired ablation site. The cylindrical tip electrode may have either straight or concave walls. Tip deflection is controlled at the proximal end by a

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INDICATIONS FOR USE

The Diagnostic/Ablation Deflectable Tip Catheter is indicated for cardiac electrophysiological mapping (stimulation and recording) and, when used in conjunction with a radiofrequency generator, for cardiac ablation.

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INSTRUCTION FOR USE

- Remove catheter from package and place in a sterile work area.

- Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.

- Connect the catheter to the recording equipment and/or the radiofrequency generator using the appropriate interface cables.

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WARNINGS

- Implantable pacemakers and Automatic Implantable Cardioverter/Defibrillator (AICDs) may be adversely affected by radiofrequency current. It is important to: have temporary external sources of pacing available during ablation; temporarily reprogram the pacing system to minimum output or OFF mode to minimize risk of inappropriate pacing. Exercise extreme caution during ablation when in close proximity to atrial or ventricular permanent leads; program the AICD to the OFF mode during the ablation procedure; and, perform complete implantable device analysis on all patients after ablation.
- Patients undergoing septal accessory pathway ablation are at risk for complete AV block which requires the implantation of a permanent pacemaker.
- During the trans-aortic approach, adequate fluoroscopic visualization is necessary to avoid placement of the catheter in the coronary vasculature. Intracoronary placement of the ablation catheter or RF

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energy application or both have been associated with myocardial infarction.

- Catheter ablation procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging. Catheter ablation should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure. Careful consideration must therefore be given for the use of the device in pregnant women.

- Do not expose catheter to organic solvents such as alcohol.

- Do not autoclave the catheter.

- Do not immerse proximal handle or cable connector in fluids; electrical performance could be affected.

- Do not scrub or twist the distal tip electrode during cleaning.

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- To avoid thromboemboli, intravenous heparin is used when entering the left heart during ablation, and many physicians prescribe aspirin, less often warfarin, for about 3 months afterward. No consensus yet exists about the need for short-term anticoagulation after ablation.
- Cardiac ablation procedures should be performed by appropriately trained personnel in a fully-equipped electrophysiology laboratory.
- Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade.
- Catheter advancement should be done under fluoroscopic guidance.
- Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- Always pull the thumbknob back to straighten the catheter tip before insertion or withdrawal of the catheter.
- When radiofrequency current is interrupted for either a temperature or an impedance rise (the set limit is exceeded), the catheter should be removed, and the tip cleaned of coagulum. When cleaning the tip electrode, be careful not to twist the tip electrode with respect to the catheter shaft; twisting may damage the tip electrode bond and loosen the tip electrode.
- Apparent low power output, high impedance reading or failure of the equipment to

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function correctly at normal settings may indicate faulty application of the dispersive electrode(s) or failure of an electrical lead. Do not increase power before checking for obvious defects or misapplication.

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- The Diagnostic/Ablation Deflectable Tip Catheter is intended for use with compatible* radiofrequency generators and accessories only.

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- Store in a cool, dry, dark place.

- The Diagnostic/Ablation Deflectable Tip Catheter is intended for single patient use only.

- Do not resterilize and reuse.

- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may

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lead to injury, illness, or death of the patient.

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- Electromagnetic interference (EMI) produced by the Diagnostic/Ablation Deflectable Tip Catheter used in conjunction with a Radiofrequency generator during normal operation may adversely affect the performance of other equipment.

- Electrodes and probes of monitoring and stimulating devices can provide paths for high frequency current. The risk of burns can be reduced but not eliminated by placing the electrodes and probes as far away as

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possible from the ablation site and the dispersive electrode. Protective impedances may reduce the risk of burns, and permit continuous monitoring of the electrocardiogram during energy delivery.

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*Laboratory compatibility test data on file

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NOTE: For catheters and cables with tip pin terminations, the pins are marked as follows: a) "D" denotes the tip (ablating) electrode, b) "THR" denotes a thermistor temperature sensor, and c) "TC" denotes a thermocouple temperature sensor. The catheter tip electrode must be switched from the electrogram recording equipment

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(via the generator controls or an external switch) to the radiofrequency generator power output for ablation (see applicable radiofrequency generator manual for instructions). To complete the electrical circuit, a dispersive pad must be connected to the reference electrode input on the generator. Circuit impedance prior to radiofrequency ablation should be approximately 100 ohms. For the temperature sensing catheter, verify that the generator displays a temperature near 37°C after the catheter is inserted in the patient and before applying RF power.

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DESCRIPTION

The EZ STEER Bi-Directional Electrodes Injection (further – Catheters) Deflectable Tip are multi-electrode devices that can be used for electrophysiological mapping and radiofrequency (RF) ablations of the heart. Standard features of this catheter include a braided bi-directional deflectable tip section, an array of four platinum electrodes which includes a 4mm tip dome, a high torque shaft, and a choice of either thermistor or thermocouple temperature sensing capabilities. Additionally, a variety of curve types are available in symmetric or asymmetric combinations, providing two 180° opposed, single planed curves. When used in conjunction with standard recording equipment, a compatible RF generator and a dispersive pad reference electrode all four of the electrodes are capable of intra-cardiac recording and stimulation while only the tip electrode can deliver RF energy.

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Currently, the available curves for both EZ STEER Bi-Directional Catheters include DD, FF, JJ, DF, and FJ. A Rocker Lever is used in the same manner for both models to deflect the tip.

CONSULT THE LOCAL DISTRIBUTOR FOR THE APPROPRIATE INTER-FACE CABLES.

INDICATIONS FOR USE

The EZ STEER Bi-Directional Catheters Deflectable Tip are indicated for cardiac electrophysiological mapping (stimulation and recording) and, when used in conjunction with a radiofrequency generator, for cardiac ablation.

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mapping and ablation procedures. All mapping and ablation procedures must be performed in a fully equipped electrophysiology laboratory.

Detailed Device Description

The EZ STEER Bi-Directional Catheter is available in a 7 French shaft size with a single 4 mm platinum tip-electrode and three platinum ring electrodes. Each of the ring electrodes is capable of performing both intracardiac recording and stimulation. The shaft and tip sections of the EZ STEER Bi-Directional Catheter are constructed from a polyurethane blended material with internal stainless-steel braiding to provide increased torsional control and lateral stability. All standard versions of the catheter have a combined usable shaft insertion length of 115 ± 5 cm.

1. Remove catheter from package and place in a sterile work area.

2. Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.

3. Connect the catheter to the recording equipment and/or the radiofrequency generator using the appropriate interface cables.

4. Advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning.

5. Use the Rocker Lever to deflect the catheter tip (Figure 1). When the lever is pulled back from neutral, the tip will deflect relative to the direction of rotation. The amount deflection is relative to the amount of lever rotation. When the lever is pushed forward the tip will deflect in the opposite direction. To straighten the tip, return the Rocker Lever to neutral position.

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


Figure 1

The handle has an adjustable friction control that allows the operator to use the Rocker Lever and deflecting tip in a “free” state or adjust the friction to where the Rocker Lever and tip curve are “locked” in place (Figure 2). This knob is located on the opposite side of the Rocker Lever. Out of the package, the knob will be in the “off” position which allows the freest movement for the lever and deflecting tip. The amount of friction increases as the Friction Control Knob is rotated clockwise until it reaches the fully “on” position.

Clockwise rotation from the “off” position increases the friction within the Deflection Mechanism.

Counter-Clockwise rotation from the “on” position decreases the friction within the Deflection Mechanism.

Figure 2

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7. Radiofrequency current may be re-applied to the same or alternate sites using the same

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catheter. However, in the event of a generator cutoff (impedance or temperature), the catheter must be withdrawn and the tip electrode cleaned of coagulum before radiofrequency current is re-applied. When the catheter is withdrawn, the tip electrode should be cleaned prior to re-insertion. A sterile gauze pad dampened with sterile saline may be used to gently wipe the tip section clean; do not scrub or twist the tip electrode or damage to the tip electrode bond may occur and loosen the tip electrode.

CURVE TYPES

Currently, the available curves for both EZ STEER Bi-Directional catheters include DD, FF, JJ, DF, and FJ (Figure 3).

Figure 3

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The EZ STEER Bi-Directional Catheter Deflectable Tip should be used only with an RF generator that has been shown to be safe and effective for cardiac ablation. Refer to the RF generator manual for detailed generator operating instructions for RF catheter ablation. If radiofrequency current is interrupted for either a temperature or an impedance rise (a preset limit is exceeded), the catheter should be removed, and the tip cleaned of coagulum. When cleaning the tip electrode, be careful not to twist the tip electrode with respect to the catheter shaft; twisting may damage the tip electrode bond and loosen

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the tip electrode. Use the appropriate Biosense Webster accessory cable to connect the EZ STEER Bi-Directional Catheter Deflectable Tip to a compatible RF generator.

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WARNINGS

- Significant x-ray exposure, can result in acute radiation injury as well as dose-related risk for somatic and genetic effects (dose = duration of the fluoroscopic imaging multiplied by x-ray beam intensity). Take all appropriate measures to minimize x-ray exposure to both patients and clinical staff.
- Pregnancy -- Careful consideration should be given to the use of this device in pregnant women.
- Ablation from within a coronary artery can cause myocardial injury and death. Adequate fluoroscopic visualization is necessary during the transaortic approach to avoid placement of the ablation catheter in the coronary vasculature.
- Stroke or myocardial infarction may occur in patients undergoing left-sided ablation procedures. Patients should be

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PRECAUTIONS

- Cardiac ablation procedures should be performed only by appropriately trained personnel in a fully-equipped electrophysiology laboratory.
- Do not attempt to operate the EZ STEER Bi-Directional Catheter or the RF generator prior to completely reading and understanding the applicable directions for use.
- The long-term risks of protracted fluoroscopy have not been established. Careful consideration must therefore be given for the use of the device in prepubescent children.
- The long-term risks of lesions created by RF ablation have not been established.
In particular, any long-term effects of lesions in proximity to the specialized conduction system or coronary vasculature are unknown.
- It is recommended that no more than seventy (70), one minute applications of energy be performed per catheter.

Catheter Compatibility

- The EZ STEER Bi-Directional Catheters are intended for use with a compatible RF generator (see attached supplement pages) and Biosense Webster accessories only.
- Read and follow the dispersive electrode manufacturer's instructions for use; the use of dispersive electrodes,

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**DISCLAIMER OF WARRANTY
AND LIMITATION OF LIABILITY
THERE IS NO EXPRESS OR IM-
PLIED WARRANTY, INCLUDING
WITHOUT LIMITATION ANY IM-
PLIED WARRANTY OF MER-
CHANTABILITY OR FITNESS FOR
A PARTICULAR PURPOSE, ON
THE PRODUCT(S) DESCRIBED
HEREIN. UNDER NO CIRCUM-
STANCES SHALL BIOSENSE WEB-
STER, INC., OR ITS AFFILIATED
COMPANIES, BE LIABLE FOR
ANY SPECIAL, DIRECT, INCI-
DENTAL, CONSEQUENTIAL, OR
OTHER DAMAGES OTHER THAN
AS EXPRESSLY PROVIDED BY
SPECIFIC LAW.**

**WITHOUT LIMITING THE FORE-
GOING, BIOSENSE WEBSTER,
INC. OR ITS AFFILIATED COMPA-
NIES, SHALL NOT BE LIABLE FOR
ANY SPECIAL, DIRECT, INCI-
DENTAL, CONSEQUENTIAL, OR
OTHER DAMAGES, ARISING OUT
OF THE REUSE OF ANY PROD-
UCT(S) LABELED FOR SINGLE
USE OR WHERE REUSE IS PRO-
HIBITED BY APPLICABLE LAW.**

Descriptions and specifications ap-
pearing in Biosense Webster, Inc.
printed matter, including this publica-
tion, are informational only and meant
solely to generally describe the product
at the time of manufacture and are not
made or given as a warranty of the pre-
scribed product in any way.

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Medical device is stored in a cool dry place at the temperature from +5 to +25 °C and air relative humidity from 20% to 80% (without condensation).

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STERILIZATION

Sterile. Sterilized with ethylene oxide gas

This device is packaged and sterilized for single use only. Do not reuse, reprocess or re-sterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or caused patient infection or crossinfection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient

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The products have come into contact with blood and / or body fluids, disposal and /

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or destruction according to SanPiN 2.1.7.2790-10 (B class).

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DESCRIPTION

- STERILE. Sterilized using ethylene oxide.
- For single use only. Do not resterilize.
- Do not use if the package is open or damaged.

The MYOSTAR Injection Catheter is a multi-electrode, percutaneous catheter with a deflectable tip and injection needle designed to inject agents into the myocardium. The tip of the Injection Catheter is equipped with a Biosense Webster location sensor and a retractable, hollow 27-gauge needle for fluid delivery. Tip deflection is controlled at the proximal end by a tubular handpiece in which a piston slides. When the piston is pushed forward with the thumbknob, the tip is deflected (curved). When the piston is pulled back, the tip straightens. The high-torque shaft allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site. Needle movement is controlled by moving a piston on a second handle, the injection handle, located proximal to the deflection handle. Push and hold the thumbknob forward on the injection handle, to retract the needle, release the thumbknob. This handle has a standard Luer fitting for connection to a syringe. The catheter interfaces with the Biosense Webster NOGASystem, which is standard recording equipment. The dead space of the catheter is approximately 0.1 cc.

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CONSULT THE LOCAL DISTRIBUTOR
OR MANUFACTURER FOR THE
APPROPRIATE INTERFACE CABLES.

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INDICATIONS FOR USE

The MYOSTAR Injection Catheter, along with related accessory devices, is intended for targeted delivery of diagnostic and therapeutic agents indicated for the myocardium. The MYOSTAR Injection Catheter should only be used under the guidance of a NOGA Cardiac Navigation System.

Injection without the cardiac navigation system is contraindicated.

All agents must be used in accordance with the specific agent manufacturer's Instructions for Use.

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CONTRAINDICATIONS

- Use of the catheter may not be appropriate for patients with prosthetic valves.
- Do not use this device in patients with active systemic infection.
- Use of the catheter is contraindicated for patients with left ventricular thrombus.

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ADVERSE REACTIONS

Procedures requiring the introduction of a catheter into the heart should not be attempted by physicians unfamiliar with the possible complications or adverse events.

A number of serious adverse reactions have been documented for intracardiac catheter procedures:

- Pulmonary embolism
- Myocardial infarction
- Stroke
- Cardiac tamponade
- Death

Possible Complications include, but are not limited to:

- Vascular bleeding/local hematomas
- Thrombosis
- Arteriovenous fistula
- Pseudoaneurysm formation
- Thromboembolism
- Vasovagal reactions
- Cardiac perforation
- Cardiac tamponade

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

INSTRUCTIONS FOR USE

- Remove catheter from package and place in a sterile work area.
- Connect the injection catheter to the Biosense Webster NOGA Cardiac Navigation System junction box using the appropriate cables.
- Prior to beginning the procedure, measure the needle as shown in Figure 1 to ensure proper needle extension.

1. Remove Aortic Arch Simulator (AAS) & Ruler packaged with the MYOSTAR Injection Catheter and place on a sterile field.

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2. Place the extended injection catheter on a flat sterile surface.

3. Insert the catheter into the AAS.
4. Push the thumbknob on the injection handle forward to extend the needle. Hold the thumbknob and measure the needle length in millimeters.
5. With the catheter in the AAS, flex the tip to 90 degrees. Measure the needle length. The needle length should measure no more than 50% of the target tissue depth.

6. With the catheter in the AAS, measure the needle length with the catheter tip straight. The needle length should measure no more than 50% of the target tissue depth.

7. Use the thumbscrew on the injection handle thumbknob to make adjustments to the needle length.

8. A counterclockwise rotation will lengthen the extended needle. Adjustments to the needle length should be made with the catheter tip straight and the needle extended. The catheter remains in the AAS during the entire needle measurement and adjustment process.

9. Remove the AAS from the MYOSTAR Injection Catheter. Place on sterile field.

10. Document the needle length in two positions while on a flat sterile surface: 1) while straight and 2) while deflected at 90 degrees (see Figure 1).

- Connect the agent-filled 1 cc Luer hub syringe to the injection port.
- With the needle extended, flush the catheter with the agent to be injected.
- Release the thumbknob on the injector handle to retract the needle.
- Inspect the distal tip of the catheter to insure that the distal needle tip is fully retracted.

WARNING: The distal needle tip should not protrude from the distal tip of the catheter. Extension of the needle during

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- Insert the catheter.

WARNING: Always pull the thumbknob back to straighten the catheter tip prior to insertion of the catheter. Failure to do so may result in tissue damage.

- Advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and the NOGA Cardiac Navigation System to aid in proper positioning.

WARNING: To avoid thromboemboli, intravenous heparin is used when entering the left side of the heart during the injection catheterization procedure.

- The catheter tip can be deflected to facilitate positioning by using the thumbknob to vary tip curvature. Pushing the thumbknob forward causes the catheter tip to bend; when the knob is pulled back, the tip straightens.

- Prior to needle extension, determine that the following criteria have been met:

1. Target tissue is not near the mitral valve annulus or apex of the heart.

2. Catheter tip is perpendicular to the myocardial wall using NOGA Cardiac Navigation System mapping and fluoroscopy.

3. Loop Stability is < 3 mm.

4. Unipolar voltage of ≥ 7 mV is observed if the target tissue is viable myocardium or < 7 mV for infarcted myocardium.

- When it has been determined that these criteria have been met, extend the needle from the catheter distal tip by moving the piston on the injection handle and holding it in that position. After extending the needle into the myocardium, ensure that:

1. Premature ventricular contraction(s) has been observed.

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2. No catheter pushback has been observed under fluoroscopy.

3. The recommended dosage is delivered via a slow and steady injection by actuating the syringe.

WARNING: Do not inject near the mitral valve annulus or apex of the heart.

Injecting in these areas may result in effusion or perforation.

CAUTION: The catheter tip should be perpendicular to the application site prior to injecting into the myocardium.

- Before attempting to reposition or maneuver the catheter, retract the needle into the distal catheter tip, by releasing the thumbknob on the injector handle.

WARNING: Always retract the needle before attempting to reposition the catheter. Failure to do so may result in tissue damage.

- AAS & Ruler are for use with a single catheter only. Do not re-sterilize or re-use MYOSTAR Injection Catheter, AAS, or Ruler.
 - Additional injections may be delivered to other sites using the same catheter.
 - Injection catheter must be flushed with agent to be injected prior to use in the patient.
 - Prior to removal of the catheter, for additional safety, confirm that the thumbknob has been pulled back completely and the needle has been retracted. Remove the catheter and dispose of it in an appropriate manner.
- WARNING:** Always pull the thumbknob back to straighten the catheter tip and retract the needle before catheter withdrawal. Failure to do so may result in tissue damage.

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If any of the following events / symptoms occur, the procedure should be temporarily halted and the patient should be reevaluated for suitability to continue with the treatment under investigation:

- Persistent complaints of chest pain;
- Complaints of cardiac pain associated with injections;
- Persistent hypotension;
- Complaints of shortness of breath;
- ICD shocks to stop ventricular tachycardia (VT);
- DC cardioversion or defibrillation for VT;
- There is any question as to the location of the catheter tip in relation to vasculature or the left ventricle (LV);
- Any unanticipated change in level of consciousness or neurological status.

The procedure should be aborted in the event that:

- Cardiac tamponade is strongly suspected or confirmed;
- Patient experiences recurrent, hemodynamically unstable sustained runs of VT > 150 beats / minute
- The patient experiences ventricular fibrillation (VF);
- Identification of thrombus in the LV or the aorta that was not previously present on echo or left ventriculogram;
- An aortic dissection is suspected or confirmed.

After the injections are completed, remove the MYOSTAR Injection Catheter at the conclusion of the endomyocardial injection portion of the procedure.

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WARNINGS

Cardiac catheterization procedures present the potential for significant X-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the X-ray beam intensity and duration of the fluoroscopic imaging. Catheterization should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure. Careful consideration must therefore be given for the use of the device in pregnant women.

The MYOSTAR Injection Catheter should only be positioned using the guidance of a NOGA Cardiac Navigation System. Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Do not use excessive force to advance or withdraw the catheter when resistance is encountered. Animal studies have shown that liquids injected into the myocardium using the MYOSTAR Injection Catheter may be injected more deeply into tissue than anticipated due to the forward jet of injected material. The injection depth may vary for various injected materials, as depth is dependent on the nature of the injected material. Viscosity of the injected material, needle extension length, and cardiac wall thickness should be considered prior to injection.

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PRECAUTIONS

- The MYOSTAR Injection Catheter should only be used by physicians trained in interventional procedures and cardiac mapping.
- Do not attempt to operate the MYO-STAR Injection Catheter, or the Biosense Webster NOGA Cardiac Navigation System prior to completely reading and understanding the applicable Instructions for Use.
- All agents delivered must be used in accordance with the specific agent manufacturer's Instructions for Use.
- Verify compatibility of the material to be injected prior to first clinical use with the MYOSTAR Injection Catheter.
- MYOSTAR Injection Catheter procedures should be performed by appropriately trained personnel in a fully equipped laboratory.
- Do not use the injection catheter for mapping. NOGA mapping should be performed with the NOGASTAR Catheter. The NOGASTAR Catheter shall be used according to its Instructions for Use.
- The MYOSTAR Injection catheter must be flushed with agent to be injected prior to use in the patient.
- Store in a cool, dry, dark place.
- The sterile packaging and catheter should be inspected prior to use. Do not use if package is open or damaged.

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prior to initiation of an injection procedure.

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