COPY CERTIFICATION BY DOCUMENT CUSTODIAN

I, Nishitaben Patel, of legal age, being first duly sworn, hereby
Name of document custodian

Swear (or affirm) that the attached reproduction of

IFB For Pipeline Vantage
Description of documents

Is a true and exact copy of the correct and complete original document.

N.s. Patel
Signature of document custodian (affiant)


A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

State of California

County of Los Angeles

Subscribed and sworn to (or affirmed) before me on this 22nd day of April,
2022 by Nishitaben Sureshb Patel, proved
to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

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INSTRUCTIONS FOR USE

Pipeline™ Flex Embolization Device Shield Technology™

70868-001 Rev. 09/14

TABLE OF CONTENTS

Pipeline™ Flex Embolization Device with Shield Technology™

English	4
Français.	5
Deutsch.	7
Italiano	8
Español	10
Svenska.	12
Nederlands.	13
Português	15
Suomi.	16
Dansk	18
Ελληνικά	20
Čeština	21
Magyar	23
Русский.	24
Polski	26
Türkçe.	28
Norsk	29
Slovenčina	30
Română.	32
Български	33
한국어.	35
عربي.	37
Symbol Glossary	38

English Instructions for Use

EN

Pipeline™ Flex Embolization Device with Shield Technology™

CAUTION

- This device should be used only by physicians trained in percutaneous, intravascular techniques and procedures at medical facilities with the appropriate fluoroscopic equipment.
- The Pipeline™ Flex embolization device with Shield Technology™ should be used by physicians who have received appropriate training for this device.
- Carefully inspect the sterile package and the Pipeline™ Flex embolization device with Shield Technology™ prior to use to verify that neither has been damaged during shipment. Do not use kinked or damaged components.
- The Pipeline™ Flex embolization device with Shield Technology™ is not to be used after the expiration date imprinted on the product label.
- Do not reprocess or resterilize. Reprocessing and resterilization increase the risks of patient infection and may compromise device performance.

DESCRIPTION

The Pipeline™ Flex embolization device with Shield Technology™ consists of a permanent implant combined with a guide-wire-based delivery system. The Pipeline™ Flex embolization device with Shield Technology™ implant is a braided, multi-alloy, mesh under woven from platinum/tungsten and cobalt-chromium-nickel alloy wires. A photograph of Pipeline™ Flex embolization device with Shield Technology™ is shown in Figure 1a and the design of the distal delivery system is shown in Figure 1b. The woven wires of the device provide approximately 30% metal coverage of the arterial wall surface area. The implant is designed for placement in a parent vessel across the neck of an intracranial aneurysm (IA). The expanded or un-constrained diameter is 0.25 mm larger than the labeled diameter.

Shield Technology™ is a synthetic polymer surface treatment and is not derived from any animal or human sources.


The tip coil is made of platinum-tungsten alloy, the proximal bumper is platinum-iridium alloy, and the tip, distal and proximal solder joints are composed of tin-silver. The protective sleeves are designed to protect the distal portion of the braid while the Pipeline™ Flex embolization device with Shield Technology™ is advanced through the micro catheter. The proximal bumper and resheathing pad allow the user to push the Pipeline™ Flex embolization device with Shield Technology™ out of the micro catheter when the delivery system is advanced. The resheathing pad also allows the user to resheath the Pipeline™ Flex embolization device with Shield Technology™ back into the micro catheter. The resheathing marker provides the user fluoroscopic visualization for the limit of resheathing the Pipeline™ Flex embolization device with Shield Technology™.

The Pipeline™ Flex embolization device with Shield Technology™ implant is mounted on a 304 stainless steel micro guide-wire approximately 200 cm long and compressed inside an introducer sheath. The Pipeline™ Flex embolization device with Shield Technology™ is designed to be delivered only through a compatible micro catheter of 0.027 inch (0.69 mm) inside diameter at least 135 cm in length.

Figure 1a: The Pipeline™ Flex embolization device with Shield Technology™

Figure 1b: The Pipeline™ Flex embolization device with Shield Technology™

- | | | |
|----------------------|----------------------|------------------------|
| 1. Tip Coil | 5. Fluorosafe Marker | 9. Resheathing Pad |
| 2. Proximal Bumper | 6. Delivery Wire | 10. Resheathing Marker |
| 3. Introducer Sheath | 7. Distal Marker | |
| 4. Braid | 8. PTFE Sleeves | |

INDICATIONS FOR USE

The Pipeline™ Flex embolization device with Shield Technology™ is intended for endovascular embolization of cerebral aneurysms.

CONTRAINDICATIONS

- Patients who have not received antiplatelet agents prior to the procedure.
- Patients with active bacterial infection.
- Patients in whom antiplatelet and/or anticoagulation therapy is contraindicated.
- The Pipeline™ Flex embolization device with Shield Technology™ should not be used alone as sole therapy for acutely ruptured aneurysms.

POTENTIAL COMPLICATIONS

Potential complications include, but are not limited to, the following:

- Adverse reaction to antiplatelet/ anticoagulation agents or contrast media
- Ischemia
- Blindness
- Mass effect
- Coma
- Neurological deficits
- Death
- Parent Artery Stenosis
- Device fracture
- Perforation
- Device migration or misplacement
- Perforator occlusion
- Dissection of the parent artery
- Post-procedure bleeding
- Distal embolization including to a previously uninvolved territory
- Ruptured or perforated aneurysm
- Embolism
- Seizure
- Groin injury
- Stroke
- Headache
- Thromboembolism
- Hemorrhage
- Transient Ischemic Attack (TIA)
- Hematoma or hemorrhage at the puncture site
- Hydrocephalus
- Vessel occlusion
- Infection
- Vessel perforation
- Intracerebral bleeding
- Vision impairment
- Chronic pain
- AV fistula
- Pseudoaneurysm

PRECAUTIONS

- Do not use in patients in whom the angiography demonstrates the anatomy is not appropriate for endovascular treatment, due to conditions such as severe intracranial vessel tortuosity or stenosis.
- Do not attempt to re-position the Pipeline™ Flex embolization device with Shield Technology™ after deployment past the resheathing marker.
- Placement of multiple Pipeline™ Flex embolization devices with Shield Technology™ may increase the risk of ischemic complications.
- The appropriate anti-platelet and anti-coagulation therapy should be administered in accordance with standard medical practice.
- A thrombosed aneurysm may aggravate pre-existing, or cause new, symptoms of mass effect and may require medical therapy.
- Do not place the Pipeline™ Flex embolization device with Shield Technology™ in patients whom a pre-existing stent is in place in the parent artery at the target aneurysm location.
- Use of implants with labeled diameter larger than the parent vessel diameter may result in decreased effectiveness and additional safety risks due to incomplete foreshortening resulting in an implant longer than anticipated.
- Person with known allergy to cobalt/chromium alloy (including major elements of cobalt, chromium, nickel, molybdenum or tungsten) may suffer an allergic reaction to the Pipeline™ Flex embolization device with Shield Technology™ implant.
- Person with known allergy to tin, silver, stainless steel or silicone elastomer may suffer an allergic reaction to the Pipeline™ Flex embolization device with Shield Technology™ delivery system.

COMPATIBILITY

The Pipeline™ Flex embolization device with Shield Technology™ is designed to be delivered through a compatible 0.027 (0.69 mm) inside diameter micro catheter. Unconstrained diameter of the Pipeline™ Flex embolization device with Shield Technology™ is 0.25 mm greater than the labeled diameter (on the packaging). Do not use the Pipeline™ Flex embolization device with Shield Technology™ in vessel diameters that are larger than the labeled diameter.

MAGNETIC RESONANCE IMAGING

Non-clinical testing has demonstrated that the Pipeline™ Flex embolization device with Shield Technology™ is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 3 Tesla or less.
- Spatial gradient field of 720 Gauss/cm or less.
- Maximum whole-body-averaged specific absorption rate (SAR) of 4.0 W/kg for 15 minutes of scanning

In non-clinical testing, the Pipeline™ Flex embolization device with Shield Technology™ produced a temperature rise of less than 0.6°C at a maximum whole body averaged specific absorption rate (SAR) of 4.0 W/kg assessed by calorimetry for 15 minutes of MR scanning in a 3 Tesla MR 750 GE Signa 20.0 System MR Scanner.

The Pipeline™ Flex embolization device with Shield Technology™ may create local field inhomogeneity and susceptibility artifacts which may degrade the diagnostic quality of the MRI images. Based on the non-clinical testing of the 5.0 mm device using standard views, the worst case maximum artifact was < 4 mm when subjected to 3.0 Tesla. Local field artifact from the Pipeline™ Flex embolization device with Shield Technology™ may decrease the accuracy of MR angiogram in assessing vessel luminal patency.

MR image quality may be compromised if the area is in the exact same area or relatively close to the position of the Pipeline™ Flex embolization device with Shield Technology™. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

DIRECTIONS FOR USE

- Using standard interventional radiographic technique, place the micro catheter tip at least 20 mm past the distal edge of the aneurysm. Gently retract the micro catheter to reduce slack in the micro catheter prior to inserting Pipeline™ Flex embolization device with Shield Technology™.

Note: It is recommended to use a heparinized saline drip to continuously flush micro catheter during Pipeline™ Flex embolization device with Shield Technology™ use.

- Choose a Pipeline™ Flex embolization device with Shield Technology™ with labeled diameter that approximates the target vessel diameter.
 - Select an appropriately sized Pipeline™ Flex embolization device with Shield Technology™ such that its fully expanded diameter is equivalent to that of the largest target vessel. An incorrectly sized Pipeline™ Flex embolization device with Shield Technology™ may result in inadequate device placement, incomplete opening, or migration.
 - The Pipeline™ Flex embolization device with Shield Technology™ foreshortens substantially (50-60%) during deployment. Take device foreshortening into account when deploying the Pipeline™ Flex embolization device with Shield Technology™.
- Choose a Pipeline™ Flex embolization device with Shield Technology™ with labeled length that is at least 6 mm longer than the aneurysm neck.
- Remove packaging hoop from the pouch and pull the distal end of the introducer sheath from the blue clip on the packaging hoop.
- Carefully remove system from packaging hoop until the delivery wire is exposed.
- Partially insert introducer sheath into the rotating hemostatic valve (RHV) at the catheter hub and close the RHV. Confirm back flush of the heparinized saline at the proximal end of the introducer sheath prior to advancing the Pipeline™ Flex embolization device with Shield Technology™ into the micro catheter.
- Advance introducer sheath into the RHV; visually confirm the tip of the sheath is seated deeply in the hub of the micro catheter.
- Secure introducer sheath to the hub by locking down the RHV tightly.
- Advance the proximal end of delivery wire until it aligns with proximal end of introducer sheath.
- Remove the introducer sheath.

Note: The delivery wire has a fluorosafe marker no further than 125 cm from the distal end.

Caution: The fluorosafe marker is only compatible with micro catheters with a minimum length of 135 cm.

- Advance the Pipeline™ Flex embolization device with Shield Technology™ into the micro catheter by pushing the delivery wire until the tip of the delivery wire aligns with the tip of the micro catheter.

Caution: If high forces or excessive friction is encountered during delivery, discontinue delivery of the device and identify the cause of the resistance, remove device and micro catheter simultaneously. Advancement of the Pipeline™ Flex embolization device with Shield Technology™ against resistance may result in damage or patient injury.

Caution: The presence of other indwelling endovascular stents may interfere with proper deployment and function of the Pipeline™ Flex embolization device with Shield Technology™.

- Once the tip of delivery system and micro catheter are aligned, verify that the Pipeline™ Flex embolization device with Shield Technology™ is in the desired location. Distal end of Pipeline™ Flex embolization device with Shield Technology™ should be placed at least 3 mm past the distal edge of the aneurysm.
- Begin to deliver the Pipeline™ Flex embolization device with Shield Technology™ using a combination of unsheathing the Pipeline™ Flex embolization device with Shield Technology™ and pushing the delivery wire simultaneously.

WARNING

Pushing delivery wire without retracting the micro catheter at the same time will cause the open end braid to move distally in the vessel. This may cause damage to the braid or vessel.

- After the distal end of the Pipeline™ Flex embolization device with Shield Technology™ has successfully expanded, deploy the remainder of Pipeline™ Flex embolization device with Shield Technology™ by pushing the delivery wire and/or unsheathing the Pipeline™ Flex embolization device with Shield Technology™. Resheathing and/or manipulation of the micro catheter, by locking down the delivery wire and moving both as a system, may facilitate expansion of the Pipeline™ Flex embolization device with Shield Technology™.

Caution: Under fluoroscopy, carefully monitor the tip coil during Pipeline™ Flex embolization device with Shield Technology™ deployment.

Caution: If the delivery wire cannot be retracted into the micro catheter, carefully remove the delivery wire and micro catheter simultaneously.

- Resheathing Instructions: During deployment of the Pipeline™ Flex embolization device with Shield Technology™, resheathing can be performed by advancing the micro catheter while pulling delivery wire.


- The Pipeline™ Flex embolization device with Shield Technology™ can be resheathed until the resheathing marker has reached the distal marker of the micro catheter (see Figure 2 below).

- The Pipeline™ Flex embolization device with Shield Technology™ is fully resheathed when the distal marker is retracted completely inside the micro catheter. The system is designed to allow for 2 full cycles of resheathing of the Pipeline™ Flex embolization device with Shield Technology™.

WARNING

Resheathing the Pipeline™ Flex embolization device with Shield Technology™ more than 2 full cycles may cause damage to the distal or proximal ends of the braid.

Figure 2. Pipeline™ Flex embolization device with Shield Technology™ (Resheathing schematic as seen under fluoroscopy, image not to scale).

1. Proximal End of the device
2. Micro Catheter
3. Micro Catheter Distal Marker
4. Resheathing Marker
5. Proximal Bumper
6. Delivery Wire
7. Resheathing Limit
8. Device Detached

16. After the entire Pipeline™ Flex embolization device with Shield Technology™ is deployed, advance the micro catheter through the device making sure not to dislodge the braid. When the micro catheter tip is distal to the Pipeline™ Flex embolization device with Shield Technology™, retract the delivery wire into the micro catheter tip.

Caution: If the catheter cannot be advanced through the Pipeline™ Flex embolization device with Shield Technology™, carefully remove the delivery wire through the Pipeline™ Flex embolization device with Shield Technology™ construct.

17. Carefully inspect the deployed Pipeline™ Flex embolization device with Shield Technology™ under fluoroscopy to confirm that it is completely apposed to the vessel wall and not kinked. If the device is not fully apposed or is kinked, consider using a balloon catheter, micro catheter, or guide-wire to fully open it.

DISPOSAL: The implant and/or delivery system should be disposed of or returned to the manufacturer per institutional guidelines.

QUESTIONS AND ANSWERS

Q If excessive friction is experienced during the insertion of delivery system at any time during the delivery of Pipeline™ Flex embolization device with Shield Technology™, what should I do?

A Carefully remove the entire system simultaneously (micro catheter and delivery system).

Q Can I retrieve the Pipeline™ Flex embolization device with Shield Technology™ if the distal end of the Pipeline™ Flex embolization device with Shield Technology™ has expanded at an undesirable location?

A Yes. A partially deployed Pipeline™ Flex embolization device with Shield Technology™ can be resheathed per resheathing instructions, step 15 in the Directions for Use.

Q Can I retrieve a fully deployed Pipeline™ Flex embolization device with Shield Technology™?

A Once fully deployed, the Pipeline™ Flex embolization device with Shield Technology™ cannot be removed. A second Pipeline™ Flex embolization device with Shield Technology™ can be deployed if needed.

Q Can I place a second Pipeline™ Flex embolization device with Shield Technology™ inside another Pipeline™ Flex embolization device with Shield Technology™?

A Yes. A second Pipeline™ Flex embolization device with Shield Technology™ can be placed inside another Pipeline™ Flex embolization device with Shield Technology™. After placing the first Pipeline™ Flex embolization device with Shield Technology™, advance the micro catheter over the delivery wire while keeping the delivery wire core across the Pipeline™ Flex embolization device with Shield Technology™. Position the micro catheter at the desired location and retrieve the delivery wire. Select a new appropriate Pipeline™ Flex embolization device with Shield Technology™ and deploy it as normal.

Caution: Placement of multiple Pipeline™ Flex embolization devices with Shield Technology™ may increase the risk of ischemic complications.

Q If there is a difference between the proximal and distal diameter, which Pipeline™ Flex embolization device with Shield Technology™ diameter do I choose?

A Choose a Pipeline™ Flex embolization device with Shield Technology™ that matches larger (typically proximal) vessel diameter to ensure proper anchoring.

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Medtronic

Pipeline™ Vantage Embolization Device with Shield Technology™

fr	Dispositif d'embolisation Pipeline™ Vantage avec Shield Technology™
de	Pipeline™ Vantage Embolisationsimplantat mit Shield Technology™
it	Dispositivo di embolizzazione Pipeline™ Vantage con Shield Technology™
es	Dispositivo de embolización Pipeline™ Vantage con Shield Technology™
sv	Pipeline™ Vantage emboliseringsinstrument med Shield Technology™
nl	Pipeline™ Vantage-embolisatiehulpmiddel met Shield Technology™
pt	Dispositivo de embolização Pipeline™ Vantage com Shield Technology™
fi	Shield Technology™ -käsitelty Pipeline™ Vantage -embolisaatiolaite
da	Pipeline™ Vantage-emboliseringsanordning med Shield Technology™
el	Συσκευή εμβολισμού Pipeline™ Vantage με Shield Technology™
cs	Zařízení pro embolizaci Pipeline™ Vantage s Shield Technology™
hu	Shield Technology™ alkalmazásával készült Pipeline™ Vantage embolizációs eszköz
ru	Устройство для эмболизации Pipeline™ Vantage с технологией Shield Technology™
pl	Urządzenie do embolizacji Pipeline™ Vantage wyposażone w technologię Shield Technology™
tr	Shield Technology™ özellikli Pipeline™ Vantage Embolizasyon Cihazı
no	Pipeline™ Vantage-emboliseringsenhet med Shield Technology™
k	Embolizačné zariadenie Pipeline™ Vantage s technológiou Shield Technology™
o	Dispozitiv de embolizare Pipeline™ Vantage cu Shield Technology™
bg	Устройство за емболизация Pipeline™ Vantage с технология Shield Technology™
et	Shield Technology™ tehnoloogiaga embolisatsiooniseade Pipeline™ Vantage
hr	Medicinski proizvod za embolizaciju Pipeline™ Vantage s tehnologijom Shield Technology™
lt	„Pipeline™ Vantage“ embolizavimo prietaisas su „Shield Technology™“
lv	Embolizācijas ierīce Pipeline™ Vantage ar tehnoloģiju Shield Technology™
sl	Pripomoček za embolizacijo Pipeline™ Vantage s tehnologijo Shield Technology™
sr	Pipeline™ Vantage sredstvo za embolizaciju sa Shield Technology™
pt-br	Dispositivo de embolização Pipeline™ Vantage com Shield Technology™
ko	Shield Technology™를 이용한 Pipeline™ Vantage 색전술 기구
vi	Thiết bị thuyền tắc mạch máu Pipeline™ Vantage với Shield Technology™
kk	Shield Technology™ технологиясы бар Pipeline™ Vantage эмболизациялау құрылғысы
ar	جهاز الانصمام Pipeline™ Vantage المزود به Shield Technology™

TABLE OF CONTENTS

Pipeline™ Vantage Embolization Device with Shield Technology™

English.....	3
Français.....	6
Deutsch.....	9
Italiano.....	12
Español.....	15
Svenska.....	18
Nederlands.....	21
Português.....	24
Suomi.....	27
Dansk.....	30
Ελληνικά.....	33
Česky.....	36
Magyar.....	39
Русский.....	42
Polski.....	45
Türkçe.....	48
Norsk.....	51
Slovenčina.....	54
Română.....	57
Български.....	60
Eesti.....	63
hrvatski.....	65
Lietuvių.....	68
Latviski.....	71
Slovenščina.....	74
Srpski.....	77
Português (Brasil).....	80
한국어.....	83
Tiếng Việt.....	87
Қазақша.....	90
العربية.....	96
Symbol Glossary.....	97

English Instructions for Use

en

Pipeline™ Vantage Embolization Device with Shield Technology™

CAUTION

- This device should be used only by physicians with a thorough understanding of angiography and/or percutaneous neurointerventional procedures.

DESCRIPTION

The Pipeline™ Vantage Embolization Device with Shield Technology™ consists of a permanent implant combined with a guidewire-based delivery system. The Pipeline™ Vantage Embolization Device with Shield Technology™ implant is a braided, multi-alloy, mesh cylinder woven with cobalt-chromium-nickel and platinum wires. An image of the Pipeline™ Vantage Embolization Device with Shield Technology™ implant is shown in Figure 1 and the design of the device is shown in Figure 2. The woven wires of the device provide approximately 30% metal coverage of the arterial wall surface area. The implant is designed for placement in a parent vessel across the neck of an intracranial aneurysm (IA). The expanded or unconstrained diameter is 0.25 mm larger than the labeled diameter. Shield Technology™ is a surface-modification that is not derived from any animal or human sources.

The Pipeline™ Vantage Embolization Device with Shield Technology™ implant is assembled on a guide-wire based delivery system that consists of a 304-stainless steel core wire and a 304L stainless steel hypotube. The implant is assembled over 304 stainless steel resheathing components. A Platinum-Iridium Restraint is distal to the resheathing components and is termed the Resheathing Marker. Refer to Figure 3 for the Resheathing Marker position.

The tip coil is made of platinum-tungsten alloy. The tip, distal, and proximal solder joints are a tin-silver. The ePTFE protective sleeves cover and protect the distal portion of the braid while the Pipeline™ Vantage Embolization Device with Shield Technology™ implant is advanced through the micro catheter. The Resheathing components allow the user to resheath the implant back into the micro catheter. The Resheathing Marker provides the user fluoroscopic visualization for the limit of resheathing the implant.

The Pipeline™ Vantage Embolization Device with Shield Technology™ implant is compressed inside an introducer sheath. The Pipeline™ Vantage Embolization Device with Shield Technology™ implant is designed to be delivered through a compatible micro catheter of either 0.021 inch (0.53 mm) or 0.027 inch (0.69 mm) inside diameter and minimum 135 cm in length. Refer to Table 1 for micro catheter compatibility for each device size.

Figure 1. The Pipeline™ Vantage Embolization Device with Shield Technology™


Figure 2. The Pipeline™ Vantage Embolization Device with Shield Technology™ delivery system and implant (not to scale)

- | | |
|----------------------|----------------------|
| 1. Tip Coil | 5. Fluorosafe Marker |
| 2. Proximal Bumper | 6. Delivery Wire |
| 3. Introducer Sheath | 7. Distal Marker |
| 4. Braid | 8. ePTFE Sleeves |

Labeled Diameter (mm)	Compatible catheter inner diameter	Labeled Lengths (mm)
2.50	0.021 inch (0.53 mm)	10, 12, 14, 16, 18, 20
2.75		10, 12, 14, 16, 18, 20
3.00		10, 12, 14, 16, 18, 20, 25
3.25		10, 12, 14, 16, 18, 20, 25
3.50	0.027 inch (0.69 mm)	10, 12, 14, 16, 18, 20, 25
4.00		10, 12, 14, 16, 18, 20, 25, 30, 35
4.50		10, 12, 14, 16, 18, 20, 25, 30, 35, 40
5.00		10, 12, 14, 16, 18, 20, 25, 30, 35, 40
5.50		10, 12, 14, 16, 18, 20, 25, 30, 35, 40, 45, 50
6.00		10, 12, 14, 16, 18, 20, 25, 30, 35, 40, 45, 50

DEVICE COMPATIBILITY

Micro catheter compatibility is defined on the product label:

The Pipeline™ Vantage 021 system is designed to be delivered through a compatible microcatheter of 0.021 inch (0.53 mm) inside diameter at least 135 cm in length. Compatibility testing has been performed with the Phenom 21 Catheter.

The Pipeline™ Vantage 027 system is designed to be delivered through a compatible micro catheter of 0.027 inch (0.69 mm) inside diameter at least 135 cm in length. Compatibility testing has been performed with the Phenom 27 Catheter.

INTENDED PURPOSE / INDICATIONS FOR USE

The Pipeline™ Vantage Embolization Device with Shield Technology™ is intended for endovascular embolization of cerebral aneurysms.

CONTRAINDICATIONS

- Patients with active bacterial infection.
- Patients in whom antiplatelet therapy (i.e. aspirin and clopidogrel) is contraindicated.
- Patients who have not received antiplatelet agents prior to the procedure.
- The Pipeline™ Vantage Embolization device with Shield Technology™ should not be used alone as sole therapy for acutely ruptured aneurysms.

PREPARATION FOR USE

- Choose a Pipeline™ Vantage device with a labeled diameter that is approximately equivalent to the target vessel diameter.
 - Select an appropriately sized Pipeline™ Vantage device such that its fully expanded diameter is equivalent to that of the largest target vessel diameter. An incorrectly sized Pipeline™ Vantage device may result in inadequate device placement, incomplete opening or migration.
- Choose a Pipeline™ Vantage device with labeled length that is at least 6 mm longer than the aneurysm neck.
 - Take device foreshortening into account when deploying the Pipeline™ Vantage device. The Pipeline™ Vantage device foreshortens 47 - 58% during deployment.
- Remove packaging hoop from the pouch and pull the distal end of the introducer sheath from the blue clip on the packaging hoop.
- Carefully remove device from the packaging hoop until the delivery wire is exposed.

WARNING

- Pre-deploying the distal end of the device prior to introduction into the micro catheter may cause damage to the distal end of the braid.
- Partially insert introducer sheath into the rotating hemostatic valve (RHV) at the micro catheter hub and close the RHV. Use a minimum flush pressure of 250 mmHg and confirm back flush of the saline at the proximal end of the introducer sheath prior to advancing the Pipeline™ Vantage device into the micro catheter.
 - Advance introducer sheath into the RHV; visually confirm the tip of the sheath is seated deeply in the hub of the micro catheter.

DIRECTIONS FOR USE

- Using standard interventional radiographic technique, place the micro catheter tip at least 20 mm past the distal edge of the aneurysm. Gently retract the micro catheter to reduce slack in the micro catheter prior to inserting the Pipeline™ Vantage device.

NOTE: It is recommended to use a heparinized saline drip to continuously flush micro catheter during Pipeline™ Vantage device use.

- Secure introducer sheath to the hub by locking down the RHV tightly.

CAUTION: Avoid deploying the device prior to introduction into the micro catheter.

- Advance the proximal end of the delivery wire until it aligns with the proximal end of the introducer sheath.

- Remove the introducer sheath.

NOTE: The delivery wire has a fluorosafe marker no further than 125 cm from the distal end.

CAUTION: The fluorosafe marker is only compatible with microcatheters with a minimum length of 135 cm.

- Advance the Pipeline™ Vantage device into the micro catheter by pushing the delivery wire until the tip of the delivery wire aligns with the tip of the micro catheter.

CAUTION: If high forces or excessive friction are encountered during delivery, discontinue delivery of the device and identify the cause of the resistance, remove device and micro catheter simultaneously. Advancement of the Pipeline™ Vantage device against resistance may result in device damage or patient injury.

CAUTION: The presence of other indwelling endovascular stents may interfere with proper deployment and function of the Pipeline™ Vantage device.

- Once the tip of delivery wire and micro catheter are aligned, verify that the Pipeline™ Vantage implant is in the desired location. The distal end of Pipeline™ Vantage implant should be placed at least 3 mm past the distal edge of the aneurysm neck.

- Begin to deliver the Pipeline™ Vantage implant using a combination of unsheathing the Pipeline™ Vantage implant and pushing the delivery wire simultaneously.

WARNING

- Pushing delivery wire without retracting the micro catheter at the same time will cause the open-end of the braid to move distally in the vessel. This may cause damage to the braid or vessel.
- Use in anatomy with severe tortuosity, stenosis or parent vessel narrowing may result in difficulty or inability to deploy the Pipeline™ Vantage device and can lead to damage to the Pipeline™ Vantage device and microcatheter. To mitigate potential problems as a result of increased delivery forces, reduce the load in the system by:
 - Unloading the microcatheter to the inner curves of vessel by pulling back on the system (i.e., the microcatheter and delivery wire together).
 - Continue unloading the system until advancement of the device (inside of microcatheter) is observed, while minimizing the distal tip movement to prevent loss of position.
 - Begin to re-advance the delivery wire while maintaining reduced load in the microcatheter. This process should be repeated until the device passes through tortuous area and the delivery force is decreased.

8. Resheathing Instructions:

During deployment of the Pipeline™ Vantage device resheathing can be performed by either:


- Advancing the micro catheter while pinning the delivery wire
 - Advancing the micro catheter while applying tension on the delivery wire
 - Advancing the micro catheter while gently pulling the delivery wire proximally
- During deployment, the point of no return/Resheathing limit is reached when the Resheathing marker aligns with the Distal marker of the micro catheter (see Figure 3). The Resheathing limit is the maximum length of the implant that can be deployed while maintaining the ability to fully resheath the device.
- The Pipeline™ Vantage device implant is fully resheathed when the distal marker is retracted completely inside the micro catheter. The system is designed to allow for a 2 full cycles of resheathing of the Pipeline™ Vantage device.
- After the distal end of the implant has successfully expanded, deploy the remainder of the implant by pushing and/or unsheathing the delivery wire. Manipulation of the micro catheter by locking down the delivery wire and moving both as a system may facilitate expansion of the implant.

CAUTION: Under fluoroscopy, carefully monitor the tip coil position during deployment of the Pipeline™ Vantage device.

WARNING

- Resheathing the Pipeline™ Vantage device more than 2 full cycles may cause damage to the distal or proximal ends of the braid.
- Resheathing the Pipeline™ Vantage device past the distal marker of the delivery system may cause damage to the distal end of the braid.

Figure 3. Pipeline™ Vantage Embolization Device with Shield Technology™ (Resheathing schematic as seen under fluoroscopy, image not to scale).

- | | |
|---------------------------------|----------------------|
| 1. Proximal End of device | 5. Proximal Bumper |
| 2. Micro Catheter | 6. Delivery Wire |
| 3. Micro Catheter Distal Marker | 7. Resheathing Limit |
| 4. Resheathing Marker | 8. Device Detached |

- After the entire implant is deployed, advance the micro catheter through the implant making sure not to dislodge the braid. When the micro catheter tip is distal to the implant, retract the delivery wire into the micro catheter tip.

CAUTION: Avoid advancing or retracting the Resheathing Marker within the implant without coverage of the microcatheter.

CAUTION: If the catheter cannot be advanced through the Pipeline™ Vantage implant, carefully withdraw the delivery wire through the implant.

CAUTION: If the delivery wire cannot be retracted into the micro catheter, carefully remove the delivery wire and micro catheter simultaneously as a system.

- Carefully inspect the deployed implant under fluoroscopy to confirm it is completely apposed to the vessel wall and not kinked/twisted. If the device is not fully apposed or is kinked/twisted, consider using a balloon catheter, micro catheter, or guidewire to fully open it.

CAUTION: It is not recommended to use the Pipeline™ Vantage delivery wire to influence apposition of the implant.

POTENTIAL COMPLICATIONS

Potential complications of the device and the endovascular procedure include or are synonymous with, but may not be limited to the following:

- Adverse reaction to antiplatelet/anticoagulation agents, anesthesia reactions such as pain, nausea, aspiration, or to contrast media such as burn sensation and organ damage or failure or due to radiation exposure such as alopecia, burns, skin reddening, ulcers, skin discoloration, cataracts, delayed neoplasia
- Access site complications such as edema, abscess, bleeding including retroperitoneal hemorrhage, tissues damage, hematoma, hemorrhage, and nerve damage
- Vascular complications such as vasospasm, stenosis, dissection, perforation, rupture, AV fistula formation, pseudo aneurysm, occlusion, thromboembolic complications including ischemia, occlusion, embolism (to unintended territory)
- Device malfunctions such as kink, stretching, friction, fracture, breakage, foreign body, misplacement, migration, inadequate deployment, premature deployment, non-detachment, reaction to device materials (such as hypersensitivity, hemolysis, fever, mutagenic effects, inflammation, granuloma, toxicity)
- Systemic complications such as infection, discomfort, pain, fever, shock, allergic reactions, organ damage, organ failure, hypertension, hypotension, arrhythmia, angina, myocardial infarction.
- Neurological deficits or dysfunctions including stroke, infarction, visual deficits, loss of vision, seizures, motor function, transient ischemic attack, headache, cranial neuropathy, confusion, emotional changes, coma
- Bleeding/ hemorrhagic complications.

- Visual complications include but are not limited to Amaurosis fugax/transient blindness, Blindness, Diplopia, Reduced visual acuity/field, Retinal artery occlusion, Retinal ischemia, Retinal infarction, Vision impairment including scintillations, blurred vision, eye floaters
- Decreased therapeutic response including need for target aneurysm retreatment
- Intra-cranial hemorrhage (including from aneurysm rupture), mass effect, brain edema, hydrocephalus
- Death

* Consult instructions for use for other therapy devices and medications for additional potential complication information.

WARNING

- Person with known allergy to cobalt/chromium alloy (including major elements cobalt, chromium, nickel, molybdenum) or platinum may suffer an allergic reaction to the Pipeline™ Vantage device implant.
- Person with known allergy to platinum alloy (including major elements platinum, tungsten, iridium), tin, silver, stainless steel or silicone elastomer may suffer an allergic reaction to the Pipeline™ Vantage device delivery system.
- Do not reprocess or resterilize. Reprocessing and resterilization increase the risk of patient infection and compromised device performance.
- Placement of multiple Pipeline™ Vantage devices may increase the risk of ischemic complications.
- Do not attempt to reposition the device after full deployment.

PRECAUTIONS

- Physicians should undergo appropriate training prior to using the Pipeline™ Vantage device in patients.
- The Pipeline™ Vantage device is intended for single use only. Carefully inspect the sterile package and device components prior to use. Do not use if sterile package or device components are damaged.
- Use the Pipeline™ Vantage device system prior to the "Use-By-date" printed on the package.
- Do not use the Pipeline™ Vantage device in patients in whom angiography demonstrates inappropriate anatomy for endovascular treatment, such as severe pre- or post-aneurysmal narrowing or severe intracranial vessel tortuosity.
- The appropriate anti-platelet and anti-coagulation therapy should be administered in accordance with standard medical practice.
- A thrombosing aneurysm may aggravate pre-existing, or cause new, symptoms of mass effect and may require medical therapy.
- Use of implants with labeled diameter larger than the parent vessel diameter may result in decreased effectiveness and additional safety risk due to incomplete foreshortening resulting in an implant longer than anticipated.
- Take all necessary precautions with patients in whom a pre-existing stent is in place in the parent artery at the target aneurysm location.
- Take all necessary precautions to limit X-ray radiation doses to patients and themselves by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors where possible.

HOW SUPPLIED

This device is supplied STERILE using ethylene oxide. This device is Non-pyrogenic.

STORAGE AND DISPOSAL

- This device should be stored in a dry place, away from sunlight.
- Dispose of device in accordance with hospital, administrative, and/or local government policy.

DIAGNOSTIC MAGNETIC RESONANCE (MR) IMAGING

Non-clinical testing has demonstrated that the Pipeline™ Vantage device is MR Conditional for single and overlapping stents up to 70 mm in length. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla or 3-Tesla, only.
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less.
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2-W/kg (Normal Operating Mode of Operation for the MR system).
- Maximum head SAR of 3.2 W/kg.

After 15-minutes of continuous scanning the Pipeline™ Vantage device is expected to produce a maximum temperature rise up to 4.15°C.

Artifact Information

In non-clinical testing, the image artifact caused by the Pipeline™ Vantage device extends approximately 20.2 mm from this implant when imaged using a T1-weighted spin echo pulse sequence and a 3-Tesla MR system.

Multilayer implant configuration of the Pipeline™ Vantage device does not affect its MRI compatibility, including temperature rise, torque, displacement, and artifact.

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parent artery at the target aneurysm location.

- Use of implants with labeled diameter larger than the parent vessel diameter may result in decreased effectiveness and additional safety risks due to incomplete foreshortening resulting in an implant longer than anticipated.

- Person with known allergy to cobalt/chromium alloy (including major elements of cobalt, chromium, nickel, molybdenum or tungsten) may suffer an allergic reaction to the Pipeline Flex embolization device with Shield Technology implant.

- Person with known allergy to tin, silver, stainless steel or silicone elastomer may suffer an allergic reaction to the Pipeline Flex embolization device with Shield Technology delivery system.

- Using fluoroscopy, carefully observe the spiral position of the Pipeline Flex Embolization Device with Shield Technology during installation.

- If the insertion wire cannot be re-drawn into the microcatheter, carefully withdraw the insertion wire at the same time as the microcatheter.

- Performing a full cycle of returning the Pipeline Flex Embolization Device with Shield Technology to the microcatheter more than twice may damage the proximal or distal ends of the implant.

- The wire of the insertion system is equipped with a radiopaque marker located no more than 125 cm from the distal end.

- The radiopaque marker is only compatible with microcatheters with a minimum length of 135 cm.

- The presence of other installed endovascular stents may interfere with the correct installation and proper operation of the Pipeline Flex Embolization Device with Shield Technology.

- If force or excessive resistance is felt during advancement of the Pipeline Flex Embolization device with Shield Technology, discontinue insertion, determine the cause of resistance, and remove the device with the microcatheter. Advancement of the Pipeline Flex Embolization device with Shield Technology despite resistance may result in damage to the device or injury to the patient.

- If the catheter cannot be advanced through the Pipeline Flex Embolization device with Shield Technology, carefully withdraw the Insertion System Wire through the Pipeline Flex Embolization device with Shield Technology coating Assembly.

- Installing multiple Pipeline Flex Embolization devices with Shield Technology may increase the risk of ischemic complications.

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- Pushing the insertion system wire without simultaneously pulling the microcatheter off the device will cause the open-ended implant to move distally in the vessel. This can damage the implant or vessel.

COMPATIBILITY WITH MAGNETIC RESONANCE IMAGING SYSTEMS

Nonclinical studies have shown the Pipeline Flex embolization device with Shield Technology to be conditionally MRI compatible. You can safely scan your device if the following conditions are met:

- Static magnetic field of 3 T or less.
- Spatial gradient field of 720 Gs/cm or less.
- Maximum average specific absorption rate of radiation by the human body (SAR) 4.0 W/kg during 15 minutes of scanning. In nonclinical studies, the Pipeline Flex embolization device with Shield Technology caused a temperature rise of less than 0.6 °C with a maximum average SAR of 4.0 W/kg in 15 minutes of MRI scans using an MRI system with a magnetic field strength of 3 T (750 GE Signa 20.0 System). The Pipeline Flex embolization device with Shield Technology can create local field discontinuities and interfere

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WARNINGS AND PRECAUTIONS WARNINGS

- Pre-deploying the distal end of the device prior to introduction into the micro catheter may cause damage to the distal end of the braid
- Pushing delivery wire without retracting the micro catheter at the same time will cause the opened of the braid to move distally in the vessel. This may cause damage to the braid or vessel.
- Use in anatomy with severe tortuosity, stenosis or parent vessel narrowing may result in difficulty or inability to deploy the Pipeline Vantage device and can lead to damage to the Pipeline Vantage device and microcatheter. To mitigate potential problems as a result of increased delivery forces, reduce the load in the system by:
- Unloading the microcatheter to the inner curves of vessel by pulling back on the system (i.e., the microcatheter and delivery wire together).
- Continue unloading the system until advancement of the device (inside of microcatheter) is observed, while minimizing the distal tip movement to prevent loss of position.
- Begin to re-advance the delivery wire while maintaining reduced load in the microcatheter. This process should be repeated until the device passes through tortuous area and the delivery force is decreased.
- Resheathing the Pipeline™ Vantage device more than 2 full cycles may cause damage to the distal or proximal ends of the braid.
- Resheathing the Pipeline Vantage device past the distal marker of the delivery system may cause damage to the distal end of the braid.
- Person with known allergy to cobalt/chromium alloy (including major elements cobalt, chromium, nickel, molybdenum) or platinum may suffer an allergic reaction to the Pipeline Vantage device implant.
- Person with known allergy to platinum alloy (including major elements platinum, tungsten, iridium), tin, silver, stainless steel or silicone elastomer may suffer an allergic reaction to the Pipeline Vantage device delivery system.
- Do not reprocess or resterilize. Reprocessing and resterilization increase the risk of patient infection and compromised device performance.
- Placement of multiple Pipeline Vantage devices may increase the risk of ischemic complications.
- Do not attempt to reposition the device after full deployment.

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PRECAUTIONS

- Physicians should undergo appropriate training prior to using the Pipeline Vantage device in patients.
- The Pipeline Vantage device is intended for single use only. Carefully inspect the sterile package and device components prior to use. Do not use if sterile package or device components are damaged.
- Use the Pipeline Vantage device system prior to the "Use-By-date" printed on the package.
- Do not use the Pipeline Vantage device in patients in whom angiography demonstrates inappropriate anatomy for endovascular treatment, such as severe pre- or post-aneurysmal narrowing or severe intracranial vessel tortuosity.
- The appropriate anti-platelet and anti-coagulation therapy should be administered in accordance with standard medical practice.
- A thrombosing aneurysm may aggravate pre-existing, or cause new, symptoms of mass effect and may require medical therapy.
- Use of implants with labeled diameter larger than the parent vessel diameter may result in decreased effectiveness and additional safety risk due to incomplete foreshortening resulting in an implant longer than anticipated.
- Take all necessary precautions with patients in whom a pre-existing stent is in place in the parent artery at the target aneurysm location.
- Take all necessary precautions to limit X-ray radiation doses to patients and themselves by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors where possible.

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Non-clinical testing has demonstrated that the Pipeline Vantage device is MR Conditional for single and overlapping stents up to 70 mm in length. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

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After 15-minutes of continuous scanning the Pipeline Vantage device is expected to produce a maximum temperature rise up to 4.15°C.

Artifact Information

In non-clinical testing, the image artifact caused by the Pipeline Vantage device extends approximately 20,2 mm from this implant when imaged using a T1-weighted spin echo pulse sequence and a 3-Tesla MR system.

Multilayer implant configuration of the Pipeline Vantage device does not affect its MRI compatibility, including temperature rise, torque, displacement, and artifact.

POTENTIAL COMPLICATIONS

Potential complications of the device and the endovascular procedure include or are synonymous with, but may not be limited to the

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Pipeline Flex embolization device with Shield Technology.

Caution: Under fluoroscopy, carefully monitor the tip coil during Pipeline Flex embolization device with Shield Technology deployment.

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15. Resheathing Instructions: During deployment of the Pipeline Flex embolization device with Shield Technology resheathing can be performed by advancing the micro catheter while pulling delivery wire.

- The Pipeline Flex embolization device with Shield Technology can be resheathed until the resheathing marker has reached the distal marker of the micro catheter (see Figure 2 below).

- The Pipeline Flex embolization device with Shield Technology is fully resheathed when the distal marker is retracted completely inside the micro catheter. The system is designed to allow for 2 full cycles of resheathing of the Pipeline Flex embolization device with Shield Technology.

WARNING

Resheathing the Pipeline Flex embolization device with Shield Technology more than 2 full cycles may cause damage to the distal or proximal ends of the braid.

Figure 2. Pipeline Flex embolization device with Shield Technology (Resheathing schematic as seen under fluoroscopy, image not to scale).

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

16. After the entire Pipeline Flex embolization device with Shield Technology is deployed, advance the micro catheter through the device making sure not to dislodge the braid. When the micro catheter tip is distal to the Pipeline Flex embolization device with Shield Technology, retract the delivery wire into the micro catheter tip.

Caution: If the catheter cannot be advanced through the Pipeline Flex embolization device with Shield Technology, carefully remove the delivery wire through the Pipeline Flex embolization device with Shield Technology construct.

17. Carefully inspect the deployed Pipeline Flex embolization device with Shield Technology under fluoroscopy to confirm that it is completely apposed to the vessel wall and not kinked. If the device is not fully apposed or is kinked, consider using a balloon catheter, micro catheter, or guide-wire to fully open it.

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9.2. DIRECTIONS FOR USE of Pipeline Vantage Embolization Device with Shield Technology coating

1. Using standard interventional radiographic technique, place the micro catheter tip at least 20 mm past the distal edge of the aneurysm. Gently retract the micro catheter to reduce slack in the micro catheter prior to inserting the Pipeline Vantage device.

NOTE: It is recommended to use a heparinized saline drip to continuously flush micro catheter during Pipeline Vantage device use.

2. Secure introducer sheath to the hub by locking down the RHV tightly.

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3. Advance the proximal end of the delivery wire until it aligns with the proximal end of the introducer sheath.

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

the delivery system may cause damage to the distal end of the braid
Figure 3. Pipeline Vantage Embolization Device with Shield Technology (Resheathing schematic as seen under fluoroscopy, image not to scale)

1. Proximal End of device
2. Micro Catheter
3. Micro Catheter Distal Marker
4. Resheathing Marker
5. Proximal Bumper
6. Delivery Wire
7. Resheathing Limit
8. Device Detached

10. After the entire implant is deployed, advance the micro catheter through the implant making sure not to dislodge the braid. When the micro catheter tip is distal to the implant, retract the delivery wire into the micro catheter tip.

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CAUTION: If the catheter cannot be advanced through the Pipeline Vantage implant, carefully withdraw the delivery wire through the implant.

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11. Carefully inspect the deployed implant under fluoroscopy to

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