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Operational documentation

Endodontic motor for root canal preparation in various designs
manufactured by the company
ChangZhou BoMedent Medical Technology Co., Ltd.

MEDICAL DEVICE NAME

Endodontic motor for root canal preparation in various designs:

- I. Marc Eco Endodontic Motor for root canal preparation
1. Marc Eco control unit
 2. Contra-angle handpiece
 3. Spray nozzle
 4. Power adapter

5. USB cable
6. Silicone holder
7. Wireless foot control (if necessary)
8. Operational manual

II. Endodontic motor for root canal preparation with built-in iRoot Pro apex locator

1. iRoot Pro control unit
2. Contra-angle handpiece
3. Contra-angle case – 3 pcs.
4. Spray nozzle
5. Connecting cable A, B
6. Lip hook – 3 pcs.
7. File holder
8. Power adapter
9. USB cable
10. Silicone holder
11. Operation manual
12. Wireless foot control (if necessary)

MEDICAL DEVICE MANUFACTURER AND PLACE OF MANUFACTURE

Full name: ChangZhou BoMedent Medical Technology Co., Ltd.

Brand name: ChangZhou BoMedent Medical Technology Co., Ltd

Address: No. 9 Changyang Road, West Taihu Science and Technology Industrial Park Changzhou
213000 Jiangsu China

Telephone: +86 0519-89881606

E-mail: luke.lu@bomedent.com Сайт:

www.bomedent.com

INTENDED USE/ INDICATIONS

Intended use: The endodontic motor is designed for the preparation of the root canals of the teeth, used in endodontic treatment. During the preparation of the root canal of the tooth, at the same time, it is possible to carry out an examination using the apex locator. For independent use of the apex locator, use a separate file holder.

Scope of application: Endodontic motor designed for use in dental offices of medical institutions.

Information about potential consumers of a medical device: The endodontic motor should only be used by qualified, appropriately trained medical professionals in the dental field.

Indications for treatment:

For root canal therapy and implantology.

CONTRAINDICATIONS AND WARNINGS

Contraindications:

Do not use to expand severely curved root canals;
Do not use to treat patients with hemophilia and with a pacemaker;
Use with caution when treating patients with cardiovascular disease, pregnant women and children.

Side effects:

– none.

Residual risks:

– none.

Description of the device

Medical electrical device with an internal power supply.
Protection against electric shock - Class II, type B working parts
Not suitable for use in an oxygen-rich environment.
Mode of operation - Continuous
Ingress protection: IPX0

The endodontic root canal preparation motor is a portable device powered by built-in lithium batteries and rechargeable via a USB cable. It is compact and comfortable in shape, which can provide a better view for doctors. The LCD displays parameters such as speed, torque, operating mode, position, etc. It is also possible to set and change the operating mode using the keys, as well as select factory settings and calibrations. The product also provides full rotation mode and root canal measurement mode that can be memorized, and provides users with a key functional mode for quick use.


Medical device delivery set

- I. Marc Eco Endodontic Motor for root canal preparation
1. Marc Eco control unit
 2. Contra-angle handpiece
 3. Spray nozzle
 4. Power adapter
 5. USB cable
 6. Silicone holder
 7. Wireless foot control (if necessary)
 8. Operation manual


Control
unit





Contra-angle handpiece

Spray nozzle

Endodontic motor for root canal preparation with integrated iRoot Pro apex locator

Endodontic motor works independently Скорость на дисплее Torque, this scale displays the torque in real time__ 	Apex locator works independently Meter Apex locator working length 0,5 value means that file tip is in the apical part or nearby. Figures on the scale do not mean millimeters* "4" icon lights up when file tip is outside apical area of the root. Digital screen File tip is outside the apical part of the root. 
Apex locator and endodontic motor work regardless of each other 	Apex locator controls motor rotation 

Functions and programming buttons

ON/OFF button

- Press this button to start.
- When the power is turned on, rotation is normal.
 - Press this button to start the motor, when pressed again, the motor will stop running.
- Built-in power indicator
- **Power indicator:** When the battery is fully charged, the blue indicator is on, when the battery is half-discharged, the blue and orange indicators are on at the same time; when the battery level is low, the orange indicator is on.
- **Charge indicator:** during charging, the indicator blinks orange, when the battery is fully charged, the indicator glows solid blue

Select button

- Torque, Rotation, Operation Mode and other settings.
 - 1) With the equipment turned on, briefly press this button to select the setting: torque, rotation or operating mode.

- ✓ Adjust the technical parameters with the “+/-” selection button. The sequence is as follows: “Torque - Rotation Mode”, the settings are automatically saved, press the “P” key to exit.
- 2) При длительном нажатии на кнопку выбора активируются другие установочные параметры.
- ✓ Press the “S” button to turn on “Sound - Reciprocal rotation - Apex locator - Bluetooth - Default settings - Left / Right hand - Auto Calibration”, use the “+/-” button to adjust each parameter. Press the “P” button to exit, or press the “S” button to set other parameters.

Programming button


- You can select program from P1-P10

Button +/-  

- Press this button to increase or decrease motor speed.
- Adjust other parameters as well.

OLED display

Rotation mode

Forward, Clockwise rotation. If the load is higher than the set torque limit, or when the apical foramen is reached, the motor will automatically rotate backward. When the load is removed, the motor will start to rotate in the normal direction.

Reverse rotation, counter-clockwise.

Optimum reversing torque, reciprocal rotation.

Battery symbols

Symbols indicate battery status.

The symbol will light up when charging.

: Battery is full or nearly full.

: 30-80 % charged

: Less than 30%

charged

: Battery is empty or has low voltage. Charge the battery.

: Symbol is flickering, immediately charge the battery.

NOTE:

The symbol indicates the remaining battery power. When the file is running under load, the symbol indicating the remaining battery power shows a lower value.

Sound symbols

: Maximum sound

: Medium sound

: Quiet sound


: Silent

Built-in apex locator

Apex locator controls motor rotation.

Only apex locator working, endodontic motor is off.


Only endodontic motor working, apex locator is off.

Apex locator and endodontic motor are working regardless of each other.


- Turn on the equipment. Briefly press the “S” button and select the required mode using the “+/-” buttons.


Apical scale option

- Thanks to this option, it is convenient for the user to adjust the approximate distance to the apical foramen.
 - Scale of designations from 15 to 27;
 - Depending on the value set by the user, it improves the visualization of the channel passage.
- To change the designation of the apical foramen, follow these steps:

- Set the apical line from 15 to 27

- Hold the “S” button to enter the mode of other settings, briefly press the “S” button until a scale appears indicating the working length of the apex locator. Adjust with the “+/-” buttons.


Wireless transfer

When treating molars, if the doctor needs to penetrate deeply into the patient's oral cavity, it is inconvenient to manually select the operating mode, in this case the user can connect the wireless pedal to the endomotor using a Bluetooth wireless connection.

Wireless data transfer settings

Setting the wireless data transfer of the device:

- The wireless transfer display functions are turned on and off. Long press "S" key to enter background mode, short press "S" key until Bluetooth icon appears. Press "+/-" to select "Bluetooth connection" or "Bluetooth off".
- For Bluetooth connections, the Bluetooth icon appears. When connected, a connection signal appears:

Attention:

- This wireless transmission module is one-to-one connection and does not support multiple connections, that is, one module can only connect to one device.
- If multiple devices are automatically connected at the same time, a connection error may occur. Please connect devices one by one as required..
- The wireless transmission function automatically saves the last state on startup. The function is disabled by default.

Battery power supply


The endomotor uses a rechargeable battery. It can be charged more than 300 times, depending on the conditions of use of the endomotor. "MEMORY EFFECT" does not apply to the lithium-ion battery. Recharging is possible at any time. 300 times is the minimum allowable number of times of charging, which will not affect battery life. When replacing the battery, you must use only batteries with the appropriate technical parameters indicated by the manufacturer on the product label.

Apex locator control (Optional)

We produce safe and easy-to-use medical equipment that simplifies the work of the doctor, therefore we recommend using the mode in which the apex locator controls the endodontic motor.

Apex locator control device settings:


- In startup state, press and hold "S" to enter other setting modes, short press "S" until you enter the apex locator control mode, and press "+/-" to select "yes" or "no".
- Long press "S" key or short press "P" key to exit setting, and all settings will be saved automatically.

Factory defaults

Installed programs can be reset to factory settings if any problems arise.

- Hold the "S" button to enter the secondary settings, press the "S" button again until you enter the "Factory reset" mode. Select "yes" by pressing the "+/-" buttons counting from "9" to "0", power off and restart the equipment to complete the factory reset.
- To exit the setting mode, hold down the "S" button, or briefly press the "P" button, all settings will be automatically saved.



NOTE:

When this function is enabled, all programs will disappear and reset to factory settings. Before resetting to factory settings, save all necessary program data, if necessary.

Adjusting for right/left handed

It can be used both with the right and with the left hand.

- Hold the "S" button to enter the secondary setting mode, press the "S" button again until the "Right / left hand adjustment" option appears, select the desired function with the "+/-" buttons.
- To exit the setting mode, hold down the "S" button, or briefly press the "P" button, all settings will be automatically saved.




Auto calibration

- Place a well-lubricated contra-angle handpiece into the equipment and place on a flat surface.
- Switch on the equipment. Hold the "S" button to enter the secondary setting mode, press the "S" button again until you enter the auto calibration mode. Press "+/-" buttons and select "yes" to enter auto calibration mode. For automatic calibration, the motor will start to rotate, slowly at first, then faster and faster.
- During calibration, the motor rotates with a contra-angle handpiece. Do not touch it until automatic calibration is complete.
- When the calibration is complete, the motor will stop and turn off after 30 seconds.
- To stop this mode, hold down the "S" button for a long time, or shortly press the "P" button to exit.

NOTES:

- The battery must be sufficiently charged before starting Auto Calibration. We recommend that you run this function after fully charged.
- Before starting auto calibration, first make sure the contra-angle handpiece is clean. If the tip gets dirty during use, the calibration may not be correct.
- During calibration, do not use a file or any other additional load.
- Do not shake the motor during calibration.
- Do not perform auto calibration while charging.

TECHNICAL SPECIFICATIONS

Technical characteristics of the product "Endodontic motor for the preparation of root canals of various versions"

Parameter	Value	
	iRoot Pro	Marc Eco
Dimensions ±10 mm	210 x 25 x 27 mm (including control unit with contra-angle	210 x 25x27mm (including control unit with

Parameter	Value	
	handpiece) Control unit 155 x 25 x 27 mm	Contra-angle handpiece) Control unit 155 x 25 x 27 mm
Weight of control unit With handpiece $\pm 10\%$	136 g	136 g
Power supply type	Battery power, Maximum capacity 1200 mAh/3.7 VDC	Battery power, Maximum capacity 1200 mAh/3.7 VDC
Input charger voltage	100 – 240V - 50/60 Hz 0,15 A max	100 – 240V - 50/60 Hz 0,15 A max
Output charger adapter voltage	DC5V/2A	DC5V/2A
Voltage deviation $\pm 10\%$	Max.	Max.
Frequency $\pm 10\%$	50-60 Hz	50-60 Hz
Charging time	Charge time no more than 95 min	Charge time no more than 95 minutes
Input voltage	DC3.7V	DC3.7V
Input power	Max7VA	Max7VA
Electric safety class	Class II	Class II
Type of components used	Type B	Type B
Working mode	Continuous	Continuous
Bluetooth module	Frequency range: 2400-2483.5 MHz Transmitting power: ≤ 20 dBm	Frequency range: 2400-2483.5 MHz Transmitting power: ≤ 20 dBm
Contra-angle handpiece		
Type	Type 2 - high speed pneumatic turbine handpiece	
Dimensions $\pm 10\%$	Length 67 mm	
Weight $\pm 10\%$	37 g	
Holding device	Button holding device	
Rotation speed $\pm 10\%$	140- 600 rpm	
Torque $\pm 10\%$	0.6-4.0 N · cm	
Speed gear ratio	1:1	
Test mandrel runout when rotating without load $\pm 10\%$	≤ 0.08 mm	
Force of insertion and extraction of the control mandrel	≥ 45 N	
Minimum fixation length $\pm 10\%$	15 mm	

Parameter	Value
Maximum length of the rotating tool $\pm 10\%$	35 mm
Shank type	Type 1
Shank size $\pm 10\%$	diameter 2,35 mm
Handpiece connection Design, dimensions and tolerances of the handpiece connection part (back part)	In accordance with ISO 3964:2016
Resistance to methods and means of cleaning, disinfection and sterilization	Minimum number of sterilization cycles – 250
Adhesion coefficient	1:1
Recommended light source	LED
Spray nozzle	
Dimensions $\pm 5\%$	25 x 9 mm
Weight $\pm 10\%$	1,5 g
Handpiece case	
Dimensions $\pm 5\%$	100 x 38 x 22 mm
Weight $\pm 10\%$	3,5 g
Lip hook	
Dimensions $\pm 5\%$	62 x 2 mm
Weight $\pm 10\%$	2,5 g
File holder	
Dimensions $\pm 5\%$	190 x 65 x 12 mm
Weight $\pm 10\%$	6 g
Silicone holder	
Dimensions $\pm 5\%$	Length 168,50 mm Height 41,55 mm
Weight $\pm 10\%$	115 g
Connecting cable A, B	
Length $\pm 5\%$	1000 mm/1000 mm
weight $\pm 10\%$	9,5 g/29,5 g
Power adapter	
Dimensions $\pm 10\%$	Length 70mm x Width 40mm
Weight $\pm 10\%$	67 g
Electric shock protection class	Class 2
Connector type	Micro USB
USB cable	
Length $\pm 5\%$	1000 mm
Wireless foot control	
Dimensions $\pm 5\%$	150 x 140 x 40 mm
Weight $\pm 10\%$	376,5 g

METHOD OF APPLICATION

Merc Eco Endodontic motor for root canal preparation

Battery charging

a) Plug the USB adapter cable into the charger main unit. (Fig. A, Fig. B)

b) After connecting to an external power source, the battery charge level will be displayed on the OLED screen, the scrolling LEDs will glow yellow.

c) After charging is complete, the scrolling LEDs will glow blue.

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NOTE:


- The handpiece can be charged directly with the USB cable even if there is no USB cable. But make sure the tip is dry.
- Make sure the equipment and power supply are operated in a dry and clean environment.
- Be careful when placing the endodontic motor in the charger. Do not place the handpiece in a wet or unsafe place.
- Do not use the power supply with other equipment.
- Charging time is approximately 90 minutes, the duration depends on the conditions in which the equipment is used, the freshness of the battery, the average temperature, etc.
- The equipment cannot be used while charging.
- A completely discharged battery will damage the battery. The user should charge the battery fully once a month or every two months if the handpiece will not be used for a long time.
- Protect the charger from dust, especially metal debris. Protect the charger carefully.
- Disconnect the equipment from the charger as soon as the battery is fully charged.
- Charging may not start in the following cases:

Ambient temperature lower than 0 °C or higher than 40 °C

- Battery voltage is insufficient
- Battery voltage abnormal

Contra-angle handpiece installation

It is possible to connect the contra-angle handpiece to the central unit of the equipment in 6 adjustable positions of the handpiece head. Adjust so that the handpiece seats line up with the motor center block, insert the handpiece head until you hear a click. To remove the contra-angle head, simply pull forward.

NOTE:

- Disconnect power to remove or attach the contra-angle handpiece.
- Use the contra-angle handpiece only as recommended by the manufacturer.
- Make sure the contra-angle handpiece is tightly connected to the center block of the endodontic motor.

Installing and extracting the file

Installation: insert the file to the head of handpiece and slightly turn the file until it clicks.

Eject: Press the button located on the back and eject the file.


NOTE:

- Turn off the power before installing and extracting a file.
- When inserting a file, when you hear a click, slightly tug on the file to make sure the file is firmly seated.
- Always clean the file shank before installation and use. Dirt getting into the file when inserting the file can result in loss of concentricity and poor clamping force.
- Do not exceed the rotation speed recommended by the file manufacturers.

Prepare to operation

1. Press and hold the POWER button to turn on the equipment.
2. Keep the programming button pressed until the program number suitable for the file used is displayed. When changing the torque set value or rotation mode and other operations, press the select button, select the setting item and press the "+/-" key to change.
3. Press the "+/-" button to change the speed.

Operation

Endodontic motor features


- Quick start
- User program
- Forward (auto reverse)
- Reverse rotation
- Optimal torque (reciprocal rotation)

• Rotation mode

Forward: clockwise rotation


When the load exceeds the torque limit, or when the apical foramen is reached, the micromotor will automatically start rotating backward. When the load is stopped, the motor will start to rotate normally.

Reverse rotation: counterclockwise movement.

Under normal conditions, clockwise rotation of the motor is quieter than counterclockwise rotation. Conventional rotary files are removed from the tooth canal when the motor is reversed. A beep is heard when using this function.

Optimum torque in reverse motion: reciprocal rotation.


load is lower than the set torque limit.
Clockwise and counterclockwise rotation

than the installer

moment, the motor is automatically required to exit the root canal. The direction of exit is opposite to the direction of normal rotation.

Or

This limit of spinning will start the

NOTE:

This function is specially designed for RECIPROCAL ROTATING FILES, the cutting side moves counterclockwise, clockwise - the file exit. If the load is higher than the formatted torque (torque cannot be set), the motor will continuously rotate clockwise to allow the file to exit the channel. To use next time, you just need to press the ON / OFF button.


Rotation mode setting

- Start the equipment. Press the mode select button to adjust the motor rotation mode.
- Select the motor rotation mode: adjust the forward rotation mode (auto reverse rotation) - reverse rotation - reciprocal rotation, set the values with the "+/-" button.
- Press the "P" button to exit, the settings will be automatically saved.

Speed, torque

- Start the installation with the equipment turned on
- Speed: Press the "+/-" button to adjust the speed ;
- Torque: Press the "S" button, adjust the torque with the "+/-" buttons;
- Press "P" button to exit, all settings will be automatically saved.

Shutdown

After completing the treatment, place the endodontic motor on the stand and hold the POWER button for more than one second to turn off the power.

- Automatic shutdown

If the button does not work, or approximately 3 minutes pass during no-load rotation, the power will automatically shut off to save power and prevent malfunction. However, if a load is applied while the motor is rotating, the power will not turn off.

- Last setting


The last setting before the power was turned off is remembered. The next time the power is turned on, the values are set to the same values as immediately before the power was turned off.

Endodontic motor for root canal preparation with integrated iRoot Pro apex locator

Features and installation

General operating functions of endodontic motor and apex locator


- Automatic rotation when immersed in the tooth canal
- Automatic deceleration when approaching the apical part
- Automatic backward rotation when reaching the apical part
- Automatic stop when leaving the root canal

- iRoot Pro automatically detects if the lip hook cable is inserted
- Connecting the mouthpiece cable means that the user proceeds to simultaneously prepare the canal and measure it;

Apex locator features

- When starting only apex locator mode, the user can connect the cable to start measuring the canal.

Endodontic motor features

- User program
- Forward (Auto reverse)
- Reverse rotation
- Optimum torque (reciprocal rotation)

Contra-angle handpiece installation

It is possible to connect the contra-angle to the central unit of the equipment in 6 adjustable positions of the handpiece head. Adjust so that the handpiece seats line up with the motor center block, insert the handpiece head until you hear a click. To remove the contra-angle head, simply pull forward.

NOTE:

- Disconnect power to remove or attach the contra-angle handpiece.
- Use the contra-angle handpiece only as recommended by the manufacturer.
- Make sure the contra-angle is tightly connected to the center block of the endodontic motor.

File installation and extraction

Installation: Insert the file into the contra-angle head and rotate the file slightly until it clicks.

Eject: Press the button located on the back and eject the file.

Charging

- Plug the USB adapter cable into the main unit's charger. (Fig. A, Fig. B)
- After connecting to an external power source, the battery charge level will be displayed on the OLED screen, the scrolling LEDs will glow yellow.
- After charging is complete, the scrolling LEDs will glow blue.

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- Make sure the equipment and power supply are operated in a dry and clean environment.
- Do not use the power supply with other equipment.
- Charging time is approximately 90 minutes, the duration depends on the conditions in which the equipment is used, the freshness of the battery, the average temperature, etc.
- The equipment cannot be used while charging.
- Do not allow the battery to be completely discharged, as this may damage it. If the equipment is not used for a long time, the battery should be fully charged once a month.
- Disconnect the equipment from the charger immediately after the battery is fully charged.

Control

Built-in apex locator

iRoot Pro endodontic motor has a built-in apex locator

Attention:

- Only the contra-angle handpiece in combination with the mouthpiece, approved by the manufacturer, gives an accurate result when measuring the length of the root canal of the tooth.
- Use only Niti endodontic files with a metal handle to determine the length of the root canal of the tooth.
- In some cases it is not possible to accurately determine the length of the root canal (see Chapter 9 "Troubleshooting" for more information)
- Electrical measurement of the length of the root canal takes place at minimum current. The values used in the iRoot Pro system are well below the values specified in the IEC 60601-1 Directive. However, in rare cases, patients experience an electrical shock when measuring apical length. If this happens, do not use this equipment to treat these patients.

Four working modes

- There are four operating modes for Apex Locator and Endodontic Motor.


Apex locator controls the rotation of the endodontic motor.

Only the apex locator works, the endodontic motor does not work.

Only the endodontic motor works, the apex locator does not work.

Apex locator and endodontic motor work independently of each other.

- Turn on the equipment. Briefly press the "S" button and select the required mode using the "+/-" buttons.

NOTE:

**ENDO-M
APEX-L**

In this mode, the endodontic motor and the apex locator operate simultaneously, without controlling each other. Upon reaching the apical portion, the apex locator only shows the apical length and does not control the endodontic motor to switch to counterclockwise and reverse rotation.

**ENDO-M
APEX-L**

In this mode, the endodontic motor and the apex locator operate simultaneously, the apex locator adjusts the endodontic motor. Endodontic motor automatically starts, automatically starts reverse rotation while measuring the length of the root canal by the apex locator.

Useful tips to exactly determine the canal length


- We recommend using gloves and rubber pads to isolate the teeth.
- Use a suction system and a cotton swab to remove excess fluid from the mouth.
- Avoid direct contact between the contra-angle file and the oral mucosa. Use a silicone case for the contra-angle handpiece.

Apical scale

- Due to this option, it is convenient for the user to adjust the approximate distance to the apical foramen.
- Scale of designations from 15 to 27;
- Depending on the value set by the user, it improves the visualization of the channel passage.

To change the designation of the apical foramen, follow these steps:

- Set the apical line from 15 to 27
- Hold the "S" button to enter the mode of other settings, briefly press the "S" button until a scale appears indicating the working length of the apex locator. Adjust with the "+/-" buttons.

NOTE:

- The settings are automatically saved.
- The apical point mark is available for all operating modes, including QS mode and mode in which the apex locator adjusts motor operation.

Apex locator control – Auto-start, Auto reverse, exit

We produce safe and easy-to-use medical equipment that simplifies the work of the doctor, therefore we recommend using the mode in which the apex locator controls the motor..


Connection

For complex measurement of apical length, use a contra-angle handpiece with silicone case.


- Connect the lip hook cable to the micro USB connector located at the bottom of the center unit. Insert the lip hook into the connector at the end of the cable. Make sure the connection is tight.
- Insert the mouthpiece into the patient's mouth (we recommend placing it on the opposite side of the tooth to be treated).
- Use a separate manually adjustable file clamp instead of the contra-angle handpiece for length measurements with an apex locator without an endodontic motor.
- Connect the apex locator cable to the micro USB connector located at the bottom of the

center unit. Insert the mouthpiece into one of the connectors and the file clip into the other connector. Make sure the connection is tight.

- Place the mouthpiece in the patient's mouth (on the opposite side of the tooth to be treated).

Comprehensive apical length measurement

Independent measurement of the apex with a locator

Procedure

- Press the "POWER" button to start the equipment.
- Plug in the lip hook cable and select the operating mode: the apex locator controls the endodontic motor, or the apex locator and the endodontic motor work simultaneously but independently of each other. When selecting the independent operation of the apex locator, insert the connecting cable to start measuring the canal.
* Parameters such as: speed and torque can be changed according to the user manual.
- Contour measurement

For a comprehensive measurement of apical length Connect the lip hook to the file

- 1) The file will rotate counterclockwise;

- 2) Warning - short beep;
- 3) Information on the screen: Scale indicating the length of the root canal, complete, information about the root canal: "-", Apex locator flashes continuously.

The apex locator works independently of the endodontic motor. Connect the lip hook to the file clip.

- 1) Warning - short beep;
- 2) Information on the screen: the scale indicating the length of the root canal is full, information about the root canal: "-", the word "APEX" flashes continuously.

D. Operation : After setting the above parameters, connecting the mouthpiece, measuring the contour, hook the mouthpiece on the patient's lip with your hand and the treatment can begin.

Motor rotation:

- 1) Place the file in the tooth canal, the file will start rotating at the set speed. Press the ON / OFF button to set the file in motion.
- 2) When the file reaches the apical foramen, the approximate value of the root canal will be 02, the motor will automatically slow down, the approximate value of the root canal will be 01, when it reaches the apical foramen, the motor will start rotating in the opposite direction.
- 3) 3) When rotating in the opposite direction, the scale value must go up to 05 for the motor to start rotating forward at the set speed and torque. The sequence of actions will be repeated.

Sound signals:

When the file is immersed, the sound signals sound from quiet to loud, when approaching the apical part, the sound signals become constant, when going beyond the apical part, the sound signals become more frequent.

Digital screen:


When the file is immersed, the relative designations of the root canal on the scale change, when approaching the apical part, the word "APEX" appears, the root canal scale is filled, the relative designation on the scale becomes "00"; outside the apical part, the word "APEX" flashes, the root canal scale is full, the relative value of the root canal is "-".

Rotation during various working modes.


Forward: clockwise rotation.

When the load exceeds the torque limit, or when the apical foramen is reached, the micromotor will automatically start rotating backward. When the load is removed, the motor will start to rotate normally.


Load lower than torque limit

Reverse rotation when loaded above the set torque limit

Reverse rotation with further loading and upon reaching the apical part. When the load stops, the motor stops and starts rotating clockwise.

Reverse rotation with further loading and upon reaching the apical part. When the load stops, the motor stops and starts rotating clockwise.

Reverse rotation: Counterclockwise movement.


Under normal conditions, clockwise rotation of the motor is quieter than counterclockwise rotation. Conventional rotary files are removed from the tooth canal when the motor is reversed. A beep is heard when using this function.

Counterclockwise rotation

Optimum torque in reverse motion: reciprocal rotation.

load
touch limit
counterclockwise in turn

load is higher,

load is higher,

but lower than the mouth of
the rutting moment. Bp

than installing
moment, the motor will automatically start moving backward
to exit the dental canal. The direction of exit is opposite to the
direction of normal rotation.

Rotation mode setting

- Start the equipment. Press the mode selector button to adjust the motor rotation mode.
- Select the motor rotation mode: adjust the forward rotation mode (auto reverse rotation) — reverse rotation — reciprocal rotation, set the values with the “+/-” button.
- Press the “P” button to exit, the settings will be automatically saved.

Speed,

torque

- Begin the installation with the hardware turned on
- Speed: Press the “+/-” button to adjust the speed ;
- Torque: Press the “S” button, adjust the torque with the “+/-” buttons;
- Press the “P” button to exit, all settings will be automatically saved.


Determination of the length of the tooth root using the iRoot Pro system and using X-ray.

Since the X-ray images show a two-dimensional image of a three-dimensional dental canal, therefore, in some cases, the data obtained using an X-ray does not coincide with the data obtained when determining the length of the tooth root by an electric method. This does not mean that the data obtained using the iRoot Pro system or X-ray is inaccurate.

These discrepancies are associated with morphological deviations. The actual position of the apical foramen is not always correctly displayed on X-ray images.

If there are cracks in the root of the tooth, it is impossible to accurately measure the length of the root canal, as the electric current penetrates the cracks in the root of the tooth.

Apical foramen is located towards the crown

Bending the root of a tooth may appear shorter on an x-ray than when measured with the iRoot Pro system, since the curvature of the root is measured using radiation.

SAFETY PRECAUTIONS

- When replacing the battery, it is necessary to use only batteries with the corresponding technical parameters indicated by the manufacturer on the product label.
- The product is specially designed for use with rechargeable batteries. Do not use manganese or alkaline batteries. Using and charging by mistake may result in fluid leakage or explosion.
- High external pressure may cause fluid leakage or explosion. Keep motor and battery (accumulator) away from fire, as this may cause fluid leakage or explosion.
- Sterilization with hot steam may cause fluid leakage and explosion.

- Mobile RF devices can interfere with the medical device.
- Do not use RF equipment nearby.
- If the equipment has not been used for a long time, check it before starting work.
- Do not touch the power cord with wet hands. Failure to do so may result in electric shock.
- Keep handpiece away from water and chemical solutions, which may cause short circuit, fire, etc.
- Do not disassemble or change the motor handpiece.
- If the battery (accumulator) has a liquid leak, or the motor is deformed or discolored, stop use immediately and contact your dealer.
- Do not charge the motor without a battery.
- Keep the charging interface away from conductive objects such as copper scrap.
- Use the original charger to charge the battery.
- When inserting the motor handpiece into the charger, make sure that there is a beep and the LED indicator shows charging (animation of the remaining battery charge). If charging is not specified, this function will not work and may result in burns due to heat generation or fluid leakage, please discontinue use and contact your dealer.
- The motor handpiece has an electronic circuitry (TORQUE LIMITER function) to prevent instrument breakage.
- Complete training is required in the use of this product to ensure patient safety.
- The product is used for dental treatment only by qualified personnel.
- Do not use a bent, damaged, deformed or inappropriate tool. It could cause injury due to sudden breakage or flying off during rotation.
- Do not use or leave the product in high temperature environments, such as in strong direct sunlight, in a car under the scorching sun, near a fire, or near a stove.
- Check the product for abnormal attenuation, vibration, noise and temperature (heat generation) before use. Discontinue use immediately and contact your supplier.
- Keep the motor away from dirt, this will affect the concentricity of the pitch angle and the determination of the torque.
- Disconnect power to the motor before changing the tilt angle or servicing the motor.
- Do not lubricate motor parts. Only lubricate the head and shank.
- If corrosive solutions are used during operation, please clean the motor handpiece immediately after operation.
- The product is intended for indoor use only.

- Follow the instructions for use of the tool to set the motor speed.
- The motor uses a rechargeable battery.
- Power consumption is low, when the power is off, please fully charge after each use, please charge at least 1 time per month for long-term use.
- When the battery is discharged, the motor will automatically shut off, leave it for a while, and then the motor can be turned on again, but this is unacceptable, it may damage and shorten the battery life.
- If the product is not used for a long period, remove the battery.
- If there is no power left at all, the battery may be damaged, replace the battery.
- Please read the instructions for use, maintenance and operation.

METHODS AND MEANS OF CLEANING, DISINFECTION AND STERILIZATION

Cleaning

Insert the spray nozzle into the contra-angle head (plug into the connector to connect to the control unit).

Insert the nozzle of the oil spray into the corresponding grease socket, spray the grease for

1-2 seconds until oil starts to flow out of the contra-angle head.

1. Lubricating the contra-angle handpiece

- Use a proven spray for lubrication or sterilization.
- Only lubricate the handpiece.
- Apply lubricant spray after each use and / or before each calibration and sterilization.

2. Cleaning the endodontic motor control unit

- Wipe it with a cotton cloth dampened with alcohol.

3. Cleaning the silicone holder

- Wipe it with a cotton cloth dampened with alcohol.

Attention:

- Wipe off excess oil before attaching the lubricated contra-angle handpiece to the control unit. Hold the end or tilt it to the correct position to drain by gravity. Install it after excess oil is completely drained.
- Hold the contra-angle firmly to prevent it from flying off under the pressure of the spray.
- Hold lubricant spray (optional) upright.
- When cleaning the endodontic motor control box, do not use solvents such as benzene and thinner.

Sterilization

Attention:

- The manufacturing plant does not disinfect parts of the endodontic motor.
- Do not sterilize the endodontic motor control unit.
- Do not wipe with unclean agents, or immerse in acidic water or sterilizing solutions.
- Parts of an endodontic motor that can be disinfected: lip hook, contra-angle sheath, contra-angle handpiece and file clip.

The lip hook, contra-angle case and contra-angle handpiece are non-toxic (in accordance with EN ISO 10993-1). They must be disinfected before use for each individual patient. Steam disinfection at high temperature and high pressure is recommended. Please refer to the following parameters for high temperature, high pressure steam sterilization:

a) Steam sterilization at 121 ° C (250 ° F) for 20 minutes (place in sterilization bag), high temperature sterilization should not exceed 135 ° C (275 ° F). The steam sterilizer must comply with EN 13060 standards.

b) The lip hook, contra-angle case, contra-angle handpiece and file holder can be re-disinfected.

High temperature and high pressure sterilization is only permissible for the lip hook, contra-angle case, contra-angle handpiece and file holder, and not permissible for the rest of the endodontic motor.

CHANGING THE BATTERY

The endodontic motor uses a rechargeable battery. It can be charged more than 300 times, depending on the conditions of use of the endodontic motor.


- "MEMORY EFFECT" does not apply to the lithium-ion battery. Recharging is possible at any time.
- 300 times is the minimum allowable number of times of charging, which will not affect battery life.
- Elements affect battery life:
- Have not been used for a long time. We recommend charging at least once a month.
- Do not charge the battery if the room temperature is over 30 degrees.

NOTE:

- Open only the battery compartment cover.
- Use only batteries specified by the manufacturer.
- Special battery: Li-ion battery, 16,500 cylinders, 3.7 V nominal
- Charging with this battery may cause damage
- Discontinue use immediately if you find a battery fluid leak.
- Do not use a damp cloth, alcohol, or other chemical cleaner to wipe the battery surface.
- Do not work with wet hands. Failure to do so may result in rust on the battery or cause a fire.
- If battery fluid gets into your eyes, rinse immediately with clean water and seek medical attention.

Battery replacement

1. Open the battery cover.
Slide it towards the charging connector.
2. Remove the old battery.
Remove the battery and then carefully pull the battery plug out of the center unit.
3. Insert a new battery.
4. Close the cover.
Slide the battery cover up by pressing lightly on it with your finger.

COMPLIANCE WITH THE NATIONAL STANDARDS OF THE RUSSIAN FEDERATION

GOST ISO 10993-1-2011	"Medical products. Assessment of biological actions of medical devices. Part 1: Evaluation and Research "
GOST R ISO 10993-2-2009	"Medical products. Assessment of the biological effect of medical devices. Part 2. Requirements for the treatment of animals "
GOST ISO 10993-5-2011	"Medical products. Assessment of the biological effect of medical devices. Part 5. Tests for cytotoxicity: in vitro methods "
GOST ISO 10993-10-2011	"Medical products. Assessment of the biological effect of medical devices. Part 10. Studies of irritant and sensitizing effects "

GOST ISO 10993-12-2015	"Medical products. Assessment of the biological effect of medical devices. Part 12. Preparation of samples and control samples "
GOST R 52770-2016	"Medical products. Safety requirements. Methods of sanitary-chemical and toxicological tests "
GOST 31214-2016	"Medical products. Requirements for samples and documentation submitted for toxicological, sanitary-chemical tests, tests for sterility and pyrogenicity "
GOST 31209-2003	"Containers for blood and its components. Chemical and biological safety requirements. Test methods" c. 5.3
GOST 31870-2012	"Drinking water. Determination of the content of elements by atomic spectrometry methods "
GOST ISO 7785-2-2011	"Dental handpieces. Part 2. Straight and contra-angle handpieces "
GOST R 50444-2020	Medical devices, apparatus and equipment. General Specification"
GOST R ISO 15223-1-2014	"The product is medical. Symbols used for marking on medical devices, labels and accompanying documents. Part 1. Basic requirements"
GOST R IEC 60601-1-2010	"Medical electrical equipment. Part 1. General safety requirements with regard to the basic functional characteristics "
GOST R IEC 60601-1-2-2014	"Medical electrical equipment. Part 1-2. General safety requirements with regard to the main functional characteristics. Parallel standard. Electromagnetic compatibility. Requirements and tests"
GOST ISO 14971-2011	"Medical products. Applying risk management to medical devices "
GOST R IEC 60601-1-6-2014	"Medical electrical equipment. Part 1-6. General safety requirements with regard to the main functional characteristics. Additional standard. Serviceability "
GOST R IEC 62304-2013	"Medical products. Software. Life Cycle Processes "
GOST 28195-89	"Assessment of software quality. General provisions"

GOST R ISO/IEC 9126-93	Information Technology (IT). Evaluation of software products. Quality characteristics and guidelines for their application "
GOST R ISO/IEC TO9294-93	Information Technology (IT). Software Documentation Management Guide "
GOST R IEC 62366-2013	"Medical products. Design of medical devices with regard to serviceability "

PACKAGING:

An endodontic motor for root canal preparation with accessories is placed in a cardboard box that serves as a consumer and transport package.

The packaging ensures safety and excludes the ingress of contaminants from the external environment, in an undamaged state the packaging protects against external influences.

SERVICE LIFE

The service life of the medical device "Endodontic motor for the preparation of root canals of various designs" is 4 years.

STORAGE, TRANSPORTATION AND OPERATION CONDITIONS

The medical device is transported by all types of transport in covered vehicles in accordance with the current transportation rules in the manufacturer's packaging at an ambient temperature of -10 to + 50 ° C, a relative humidity of 10 to 85% (no condensation) and an atmospheric pressure of 50 up to 106 kPa.

Transportation without the manufacturer's packaging does not guarantee its safety. Damage resulting from transportation without the manufacturer's packaging is eliminated by the consumer.

A medical device released from transport packaging must be stored at an ambient temperature of -10 to + 50 ° C, a relative humidity of 10 to 85% (non-condensing) and an atmospheric pressure of 50 to 106 kPa. Do not leave the room near fire fighting equipment and flammable materials.

A medical device must be used at an ambient temperature of +10 to + 40 ° C, a relative humidity of 10 to 70% (no condensation) and an atmospheric pressure of 70 to 106 kPa.

MAINTENANCE AND ROUTINE REPAIR

The manufacturer guarantees the high quality of the equipment under the conditions of its correct installation, operation and maintenance. The warranty does not cover batteries as they are consumable parts.

You can find detailed information in the warranty card.

ERROR CODES

	Error code	Error	Reason	Check and troubleshooting
While working	E-0	Self-diagnostics error	Circuit malfunction	Contact the dealer
	E-1	Overload	Motor handpiece is blocked (auto reverse mode)	Release the load
	E-2	Overvoltage	Circuit malfunction	Contact the dealer
While charging	E-3	Charging device is broken	Circuit malfunction	Contact the dealer
	E-4	Low battery voltage	Battery is broken	Change the battery
	E-5	High battery voltage	Battery is broken	Change the battery
	E-8	Low adapter voltage	Adapter is broken	Change the adapter
	E-9	High adapter voltage	Adapter is broken	Change the adapter
During calibration	E-7	Calibration error	Motor is damaged or handpiece is blocked	Change the motor or handpiece

TROUBLESHOOTING

If the equipment is not working properly, this does not always mean that the central unit is defective. Please review the list below first to rule out user error or other reasons before contacting the seller. If the problem persists, please contact your dealer or manufacturer.

Problem	Reason	Solution
Fail to switch on	Battery is out of charge.	Charge the battery.
	Battery is taken out.	Insert the battery.
	Inner fuse is burnt.	Contact the dealer.
No charge	The battery is fully discharged.	Replace with a new battery.
	The power plug is not plugged into the outlet.	Insert the power plug into the outlet.
	Charger is contaminated	Clean the charger
	The digital screen displays an error.	Look for the error code in Chapter 8.
	The handpiece panel does not show anything even when connected to a charger.	Contact the dealer.

	Charger is not working	Contact the dealer.
Motor handpiece fails to rotate	Contra-angle handpiece is blocked	Clean or re-insert the handpiece
Sound signal while motor is working	Contra-angle handpiece head is stuck	Clean the handpiece head
Sound signal while motor is working	Select auto reverse rotation mode, remove the contra-angle, no clicking sound when working	Contact the dealer
	Select auto reverse rotation mode, remove the contra-angle, there is a clicking sound during operation	Clean the handpiece head
	The rotating shaft of the contra-angle handpiece is dirty	Clean the handpiece head
Apex locator fails to switch on	Incorrect or insufficient contact with the measuring cable, lip hook, file holder	Reconnect the cables
	Connecting cable, file holder are worn	Replace the connecting cable, file holder
Apex locator doesn't control motor operation	Check if the working mode is set properly	Carefully read the instructions, make sure that the operating mode is set in which the apex locator regulates the operation of the endomotor
	If the measuring cable and lip hook are well connected	Read the instructions carefully, make sure to install and connect all connections properly
	If the contra-angle handpiece is in direct contact with the human body	Read the manuals carefully. Put the case on the contra-angle handpiece
	Check if the contra-angle handpiece is installed properly	Read the manual carefully, make sure the handpiece is installed properly.
Warning signal fails to appear	Check if the sound is on	Turn on the sound

DISPOSAL

The endodontic motor for root canal preparation does not contain harmful substances and components that pose a danger to human health and the environment during and after the end of its service life. After the expiration of the service life or write-off as a result of failure, they must be disposed of. Disposal should be carried out as waste of class A according to SanPiN 2.1.7.2790-10 (epidemiological safe waste, in composition close to solid waste). The lithium battery must be disposed of in accordance with applicable laws.

WARRANTY

The Endodontic Root Canal Preparation Motor (including the center unit and contra-angle handpiece) is guaranteed for 1 year

The manufacturer guarantees the compliance of the endodontic motor with the requirements established by the Operational Documentation, provided that the consumer observes the conditions of operation, transportation and storage.

The warranty period for the storage of the medical device (including the central unit and contra-angle handpiece) is 1 year (excluding the battery and consumables).

When buying a product, ask for the correct filling of the warranty card:

affixing the seller's stamp and the date of sale. The warranty period for the product starts from the date of purchase. In the absence of such a mark, the warranty period is calculated from the date of manufacture of the product.

In case of detection of faults in the product during the warranty period, the buyer can contact the seller for its repair only if there is a warranty card.

The warranty applies only to those cases when the product is out of order through no fault of the buyer!

Delivery to the service center and back is carried out at the expense of the client.

The warranty does not apply under the following circumstances:

- a) When using the product, the requirements and recommendations described in the operational documentation were not followed;
- b) The user has independently opened and disassembled the product;
- c) There are no documents for the purchase of this product.

Information about undesirable events that have signs of an adverse event (incident) is accepted by an authorized representative of the manufacturer.

"Novgodent", LLC

355003, Stavropol, Mira Street, 367/21

E-mail address of the corporate entity: accident@mail.ru

Tel.: + 7 8652-525-888

VERSION OF THE OPERATIONAL DOCUMENTATION

Version 1, 2021


MEDICAL DEVICE MARKING

Interpretation of symbols used in the marking of consumer and transport packaging - cardboard box on a medical device

"Endodontic motor for the preparation of root canals of various designs"

Manufacturer

	Authorized Representative in the European Community
	Serial number
	Date of manufacture
	The product fully complies with the requirements of the EU Directive 93/42 / EEC on medical devices
	Follow the instruction manual
	Read the documents supplied with the equipment (User's Manual)
	Type B device
	Direct current
	Proper disposal of the product
	Product of class II electrical safety
	For indoor use only

ELECTROMAGNETIC COMPATIBILITY

An endodontic motor for root canal preparation requires special measures to ensure electromagnetic compatibility and must be installed and commissioned in accordance with the EMC information in the OPERATING DOCUMENTATION;

b) an indication that the use of mobile RF communications equipment may affect MEDICAL ELECTRICAL EQUIPMENT.

Table 1 – Instruction and manufacturer's declaration – ELECTROMAGNETIC EMISSION

Guidance and Manufacturer's Declaration - Electromagnetic Emissions
--

The root canal preparation endodontic motor is intended for use in the electromagnetic environment defined below. The purchaser or user of the Endodontic Root Canal Preparation Motor should ensure that it is used in the specified electromagnetic environment.		
Electromagnetic emission test	Compliance	Electromagnetic conditions - guidance
CISPR 11 radio interference	Group 1	The [ME EQUIPMENT or ME SYSTEM] uses radio frequency energy only for its internal function. RFI emissions are low and likely will not cause any malfunction in nearby electronic equipment
CISPR 11 radio interference	Class B	Mercury® T-Fine-II Endodontic Motors are suitable for use in any location other than residential buildings and buildings directly connected to the electrical distribution network supplying residential buildings. Can be used in residential buildings and buildings directly connected to the electrical distribution network supplying residential buildings with the following warning: Warning. This equipment / system is intended for use by healthcare professionals only. This equipment / system may cause poor radio reception and disrupt the operation of nearby equipment. In this case, it may be necessary to take mitigation measures such as reorienting, relocating the Mercury® T-Fine-II Endodontic Motor.
Harmonic current components according to IEC 61000-3-2	Class B	
Voltage fluctuations and flicker according to IEC 61000-3-3	Complies	

Table 2 – Instruction and manufacturer's declaration – INTERFERENCE IMMUNITY - for Endodontic motor for root canal preparation

Instruction and manufacturer's declaration – interference immunity

The root canal preparation endodontic motor is intended for use in the electromagnetic environment defined below. The purchaser or user of the Endodontic Root Canal Preparation Motor should ensure that it is used in the specified electromagnetic environment.

Interference immunity test	Test level according to IEC 60601	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) according to IEC 61000-4-2	±6 κB – contact discharge ±8 κB – air discharge	±6 κB - contact discharge ±8 κB – air discharge	There are no special instructions for covering the floor and the room, the relative air humidity is less than 80%
Nanosecond impulse noise according to IEC 61000-4-4	±2 κB – for power lines ±1 κB – for input/output lines	±2 κB – for power lines ±1 κB – for input/output lines	Mains power quality should be consistent with a typical commercial or hospital environment.
High energy microsecond impulse noise according to IEC 61000-4-5	±1 κB when applying interference through the "wire-wire" scheme " ±2 κB when applying interference through the "wire-ground" scheme	±1 κB when applying interference through the "wire-wire" scheme " ±2 κB when applying interference through the "wire-ground" scheme	Mains power quality should be consistent with a typical commercial or hospital environment.
Voltage lacks, short interruptions and voltage changes in input power lines according to IEC 61000-4-11	<5% U_H (voltage fail >95% U_H) during 0,5 period 40% U_H (voltage fail 60% U_H) during 5 periods 70% U_H (voltage fail 30% U_H) during 25 periods <5% U_H (voltage fail >95%	<5% U_H (voltage fail >95% U_H) during 0,5 period 40% U_H (voltage fail 60% U_H) during 5 periods 70% U_H (voltage fail 30% U_H) during 25 periods <5% U_H (voltage fail >95%	The quality of the power consumed from the mains must correspond to the standard commercial or hospital environment.

	U_H) during 5 s	U_H) during 5 s	
Power frequency magnetic field (50/60 Hz) according to IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic field levels should be maintained in accordance with typical operating conditions
Note - U_H is the voltage level of the electrical network before the test action is applied.			

Table 3 – Guidance and manufacturer's declaration - IMMUNITY - for non-LIFE SUPPORT ME EQUIPMENT and ME SYSTEMS

Guidance and Manufacturer's Declaration – Interference Immunity			
The root canal preparation endodontic motor is intended for use in the electromagnetic environment defined below. The purchaser or user of the Endodontic Root Canal Preparation Motor should ensure that it is used in the specified electromagnetic environment.			
Interference immunity test	Test level according to IEC 60601	Compliance level	Electromagnetic environment - guidelines
			The distance between the mobile radiotelephone systems in use and any component of the Mercury® T-Fine-II Endodontic Motor, including cables, must be at least the recommended separation calculated using the expressions below for the frequency of the transmitter.
Conducted disturbances by radio frequency electromagnetic fields according to IEC 61000-4-6	3 (rms) from 150 kHz to 80 MHz	B $[V_1]$, B	Recommended separation distance: $d = \left[\frac{35}{V_1} \right] \sqrt{P}$

RF electromagnetic field according to IEC 61000-4-3	3 V / m from 80 MHz to 2.5 GHz	$[E_1], \text{V/m}$	$d = \left[\frac{35}{E_1} \right] \sqrt{P}$ (80-800 MHz); $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ (800 MHz - 2.5 GHz), where d is the recommended separation distance, mb; P is the rated maximum output power of the transmitter, W, established by the manufacturer. Field strengths from fixed RF transmitters, as determined by observations of the electromagnetic environment a), should be lower than the compliance level in each frequency band b). Interference can occur in the vicinity of equipment marked with the 
a) Field strengths due to propagation of radio waves from fixed radio transmitters such as base stations of radiotelephone networks (cellular / wireless) and land mobile radio stations, amateur radio stations, AM and FM broadcasting transmitters, television transmitters cannot be calculated with sufficient accuracy. For this, practical field strength measurements must be made. If the measured values at the site of the Root Canal Preparation Endodontic Motor exceed the applicable compliance levels, the operation of the Root Canal Preparation Endodontic Motor should be monitored to verify that it is functioning properly. If the observation reveals an abnormality, then it may be necessary to take additional measures, such as reorienting or moving the Endodontic Motor for root canal preparation.			

b) Outside 150 kHz to 80 MHz, field strengths should be less than, V / m. Notes (edit)

1 At 80 MHz and 800 MHz, the higher field strength applies.

2 Expressions may not apply in all cases. The propagation of electromagnetic waves is affected by absorption or reflection from structures, objects and people.

Warranty card

Dear customer:

The warranty covers:

1. The iRoot Pro / Marc Eco medical device (including control unit and contra-angle handpiece) is guaranteed for 1 year (excluding battery and accessories).
2. The hardware warranty does not apply under the following circumstances:
 - a) When using the equipment, the requirements and recommendations described in this Manual were not followed;
 - b) The user has independently opened and disassembled the equipment;
 - c) There are no documents for the purchase of this equipment.
3. Fill in the following information, then send to us along with the equipment.

Customer's name: _____

Telephone number: _____

Address:

Defect description:

(Also describe: When, where and how it happened how many times it happened)

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(Translation)

(2021) SCCZZ Zi, No. 14070

Applicant: Changzhou Bomedent Medical Technology Co., Ltd., Domicile: 3rd Floor, Building B1, No. 9 Changyang Road, West Taihu Science and Technology Industrial Park, Changzhou.

Legal representative: Li Yanlong, male, born on October 28, 1972, ID No. 340402197210281213.

Agent: Li Fei, male, born on June 11, 1990, ID No. 341222199006113575.

Issue under notarization: Seal

This is to certify that the seal of "Changzhou Bomedent Medical Technology Co., Ltd." affixed on the Operational documentation Endodontic motor for root canal preparation in various designs manufactured by the company Changzhou Bomedent Medical Technology Co., Ltd. attached hereto presented by Li Fei, the agent of Changzhou Bomedent Medical Technology Co., Ltd., is found to be authentic.

Notary: Zhou Junyi (Signature)

Changzhou Notary Public Office (Seal)

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