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
Bio-Medicus Life Support™

Catheter and Introducer - Extracorporeal Membrane Oxygenation and Extracorporeal Life Support

1. Product description

The Bio-Medicus Life Support catheter is a single-lumen catheter used to drain or reinfuse blood. The introducer facilitates proper insertion and placement of the appropriately sized catheter over a guidewire within the vessel. These devices are intended to perfuse vessels or organs in conjunction with extracorporeal support, including Extracorporeal Membrane Oxygenation (ECMO) and Extracorporeal Life Support (ECLS). Catheter models with tip lengths of 18 cm (7.09 in), 50 cm (19.7 in), or 55 cm (21.7 in) and without additional side holes may be used as either drainage or reinfusion catheters.

These devices are sterile, nonpyrogenic, disposable, intended for single use only, and have been tested for up to 21 days of use. Do not store the product above 25°C (77°F).

2. Indications for use

The Bio-Medicus Life Support catheters and introducers are single-lumen drainage or reinfusion peripheral-access catheters to be used in ECMO or ECLS with an extracorporeal circuit intended for use in adult and pediatric patients with acute respiratory failure or acute cardiopulmonary failure, where other available treatment options have failed, and continued clinical deterioration is expected or the risk of death is imminent.

3. Contraindications

Alone, the catheter and introducer are not medical treatment devices. The introducer is to be used only with the appropriately sized Bio-Medicus Life Support catheter. These devices are intended for use only as indicated in these instructions for use. Do not insert the catheter in a vessel that has arterial dissection or severe peripheral atherosclerosis.

4. Warnings and precautions

- Read all warnings, precautions, and instructions for use carefully before use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.
- Verify the use-by date.
- This device was designed for single use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Do not autoclave.
- Do not use if the sterile package has been damaged or if the seal has been broken in any way.
- Do not use alcohol or alcohol-based fluids for lubrication.
- Thoroughly inspect the catheter before use to verify that the lumen and side holes (if applicable) are patent and that the catheter has not been damaged or kinked.
- Use caution when assembling or manipulating the introducer within the catheter, or moving the suture ring.
- Avoid false lumen catheterization.
- Never advance, position, or place the catheter without the introducer. This action may kink and damage the catheter.
- Use care to protect the side hole areas from kinking as they are passed into the vessel.
- Manipulate the catheter to prevent kinks or any restrictions that may alter flow.
- Manipulation of the heart with the catheter in place may cause flow restrictions.
- Avoid kinking when attaching tubing and when anchoring or untethering the catheter body.
- Use caution when attaching the catheter to other devices in the extracorporeal circuit to prevent damage to the device being connected or the catheter.
- Clamping too close to the connector may damage the connector end.
- Do not clamp wire-reinforced sections. Clamping on reinforced sections may result in permanent wall distortion or lumen collapse.
- Since the catheter wall is opaque due to the presence of the spring, be sure to completely remove the air from the catheter during the priming procedure.
- Avoid overinsertion of the catheter and impingement of the catheter tip against the vessel, ventricular, or atrial walls.
- The use of cold solutions and cardioplegic agents, administered externally or via the coronary circulation, may increase the stiffness of urethane plastic. Increased stiffness may contribute to ventricular perforation under conditions of catheter tip pressure during manipulation of the heart.
- Do not attempt to cross-clamp the catheter with the introducer in place.
- Carefully monitor for signs of occlusion on both the inlet and outlet catheters.

Instructions for Use English 9

master:20190404, alldocMod:20180927, alldocMod2:20171221, alldocStyle_graph:20190719, alldocStyle:20180308, templates_global:20180215, master_pageset:20180222, pagespecs:20180129, body_extnt:20171025, languageLookup:20171025, table:20180606, characterJoin:20180903, 3of9:20110211, 3of9_svg:20110211

- The user assumes responsibility for any variations in extracorporeal procedures that could compromise the intended use of this device.
- Initiate adequate anticoagulation therapy before patient catheterization.
- Use of this product beyond the duration established by the available in vitro or in vivo data has not been demonstrated. Refer to Section 8 for a summary of bench and animal studies performed. Refer to Section 6 for possible clinical observations that may necessitate or predict the need for device replacement throughout the duration of use.
- Only physicians trained and experienced in using percutaneous catheterization techniques (such as the Seldinger technique), ECMO, and ECLS should use this device.
- Ischemia may occur distal to the catheter placement site.
- Catheter tip placement should be confirmed by using fluoroscopy.
- Use extreme caution to ensure that the return catheter and drainage catheter are not confused while being connected to the perfusion circuit.
- In some cases, suture rings may not be detectable on x-ray due to the variances in anatomical densities and x-ray techniques. Modifying the x-ray angle or distance may aid in detection.
- When utilizing vacuum for augmented drainage, do not exceed -100 mmHg.
- Maximum pressure is 300 mmHg.
- Do not expose to MRI while the catheter is inserted in the patient.

5. Adverse effects

The following potential adverse effects are associated with the use of the product: air embolism, anaphylaxis, death, delayed procedure, foreign-material embolism, hemolysis, hypercarbia, hypoxia, infection, irritation, ischemia, major blood loss, minor blood loss, organ dysfunction, sensitization, thromboembolism, toxicity, and vessel damage. Vessel damage and complications at the puncture site may occur if the instructions for use are not followed.

6. Complications

Possible complications include those normally associated with ECMO, ECLS, anticoagulation, and large-bore venous and arterial catheterization.

These clinical observations may necessitate or predict the need for device replacement throughout the duration of use:

- Decreased blood flow due to a clot, obstruction, or other factor that cannot be managed by catheter repositioning
- Significant or uncontrolled blood loss from the catheter
- Entrained air in the catheter or connectors
- Mechanical injury to access vessels (perforation, tearing, dissection), injury at or above the site of insertion, and tissue damage during the insertion and placement of the catheter
- Mechanical failure due to kinking of the catheter, or material failure due to incorrect suturing
- Hemorrhage related to bleeding, catheter malposition, or dislodgement around the insertion site
- Heart, vessel, or lung damage; anemia; liver or kidney failure; and stroke

Note: The benefits of catheterization for extracorporeal circulation must be weighed against the risk of systemic anticoagulation and subsequent propensity for hemorrhage.

7. Sterilization

All Bio-Medicus catheters are sterilized (EtO) before shipment. The lot number is designated on the outer package label. Do not use this device if the sterile package has been damaged or if the seal has been broken. Resterilization is not recommended under any circumstances as this device is intended for single use only.

8. Testing Summary

- An in vivo study was conducted in an ovine model for 7 days to determine the thrombotic potential of these catheters. The study was conducted with systemic heparin with ACT adjusted to 1.5 times normal. After 7 days at systemic low-flow rates (1-2 L/min), there were no visible clots at any location in these devices. Thrombus formation is a risk in all ECLS procedures. Anticoagulation in the in vivo model reflects continuous heparin drip administration, which is typical of all ECLS platforms.
- Bench studies were performed to demonstrate biocompatibility of the device materials with prolonged exposure to circulating blood. The device materials met all biocompatibility requirements.
- Bench studies were performed under simulated in vitro use conditions using blood analog at body temperature (37°C/98.6°F) to demonstrate device functionality, performance, and durability in compliance to performance specifications. Device preconditioning included extended exposure (21 days) to simulated in vivo use conditions to demonstrate safety and reliability. All devices met performance requirements.

9. Instructions for use

Note: The introducer has been designed to accommodate a guidewire (not supplied with the device):

- 0.635 mm (0.025 in) diameter guidewire (catheter sizes 9 Fr (3.0 mm), 11 Fr (3.7 mm), 13 Fr (4.3 mm), 15 Fr (5.0 mm))
- 0.97 mm (0.038 in) diameter guidewire (all other sized catheters)

1. Select a catheter size that is appropriate for the patient based on criteria including desired flow rate, body surface area, vasculature size, and catheter flow charts. Various sizes are available in order to select the appropriate catheter for each specific procedure. (Table 1 and Table 2 at the end of this document may be used to assist in catheter selection.)

Note: Flow curve data shown in the Figures section were obtained via bench testing using room-temperature water. Other factors for a given patient could result in different outcomes than the data presented.

2. Remove the catheter and introducer from the shipping carton and inspect the protective, sterile package for damage.

3. Using aseptic technique, open the protective, sterile package containing the catheter and introducer.
Note: Use a Bio-Medicus™ insertion kit when performing Seldinger or percutaneous catheterization technique.

Note: For catheter sizes 9 Fr (3.0 mm), 11 Fr (3.7 mm), 13 Fr (4.3 mm), 15 Fr (5.0 mm): select either the standard solid-core introducer for traditional placement method, or the hollow-core guidewire-capable introducer for Seldinger technique placement method. Use an insertion kit that includes a 0.635 mm (0.025 in) guidewire.

4. Insert the introducer through the hemostasis cap and ensure the introducer hub is fully inserted into the hemostasis cap.

5. Assemble the catheter and introducer. Insert the introducer into the lumen of the catheter through the connector and advance it until the hemostasis cap is pushed firmly onto the barb connector.

Note: Tighten the luer fitting cap on the catheter connector (if present) to ensure that the luer fitting is completely sealed.

6. Insert the catheter using a guidewire placement Seldinger technique. Determine the length of catheter to be introduced and note the insertion markings on the catheter.

Note: Some catheters are labeled with length markings for depth reference during catheter placement.

Note: Confirm the anatomical location of the guidewire before inserting the catheter. Anatomical location of guidewire is critical before catheter insertion.

Caution: Ensure that the catheter size selected is of adequate size for the vessel to allow distal perfusion of the limb when the catheter is in place. Improper catheter size may be difficult to advance. The vessel must be large enough to ensure perfusion and venous return.

Caution: Do not kink guidewire during insertion, vessel dilation, and manipulation or it may be difficult to remove.

7. Position the adjustable radiopaque suture ring to the desired depth.

8. Pass the catheter and introducer assembly over the guidewire, advancing until the guidewire can be grasped beyond the introducer hub.

9. Insert tip of the catheter introducer into the vessel approximately 3 cm (1.18 in) to 5 cm (1.97 in). Grasp the guidewire and hold stationary while advancing the catheter to the desired depth. Blood leakage from the proximal drainage holes may occur during catheter insertion until these holes are fully advanced into the vessel.

Note: Advancing the catheter and introducer with a rotating motion will make catheter placement easier.

Note: When the bold band closest to the introducer hub is in alignment with the edge of the hemostasis cap, the introducer tip and catheter tip are at the same depth (Refer to Figure 1). For drainage catheter sizes 15 Fr to 29 Fr (5.0 mm to 9.7 mm), the alignment of the catheter and introducer ends is indicated by the tip of the triangle-shaped symbol.

10. **For drainage catheter sizes 21 Fr to 29 Fr (7.0 mm to 9.7 mm) only:** When advancing the venous catheter into the superior vena cava (SVC), verify positioning with tactile feel, indirect visualization method, or echocardiographic imaging.

Note: For drainage models, graduated marks on the hub end of the introducer allow users to approximate the introducer length beyond the tip of the catheter.

11. For proper intracardiac advancement, once the side holes are intravascular, withdraw the introducer to the base of the triangle-shaped symbol representing the tapered tip of the introducer.

Warning: Do not withdraw to the point of the triangle-shaped symbol because of the risk of SVC, innominate vein perforation or vascular, cardiac trauma during intravascular advancement. Never force the catheter. It should advance easily.

12. When desired catheter placement has been achieved, remove the guidewire and retract the introducer clear of the clamp site.

Caution: To maintain hemostasis, do not completely remove the introducer or hemostasis cap until the clamp is applied.

Warning: Ensure that the clamp is properly positioned on the catheter tubing between the tubing connector and wire-reinforced sections. Clamping too close to the connector may damage the connector end. Clamping on wire-reinforced sections may result in permanent wall distortion or lumen collapse.

Warning: When the bold band closest to the vascular end of the introducer is visible at the edge of the hemostasis cap, the catheter can be safely cross-clamped. Do not attempt to cross-clamp the introducer.

Note: Perfusion tubing clamps must be used. Many clamps are too weak to ensure complete clamping.

13. Remove the hemostasis cap, exposing the barbed connector.

Note: Drainage catheter caps are color-coded blue.

Note: The 18 cm (7.09 in), 50 cm (19.7 in), and 55 cm (21.7 in) insertion-length catheters (with side holes at the tip only) may be used as either a return catheter or a drainage catheter. The hemostasis caps for these catheters are clear, not color-coded.

14. Slowly fill the remainder of the catheter tube with sterile priming fluid, displacing any air out through the tube connector end before connecting fluid-filled tubing.

15. Attach the circuit, taking care to connect the drainage catheter to the drainage line (inlet to pump) and the return catheter to the return line (outlet from oxygenator) to ensure an air-free connection. The tubing may be further secured with cable ties at the barbed connector.

Note: The drainage catheter is identified by an arrow in the clamp area next to the Bio-Medicus Life Support logo, indicating flow out of the connector. See Figure 1.

16. Secure the device to the patient following standard medical practice.

Note: In order to avoid accidental displacement during extracorporeal circulation, suture may be wrapped around the suture ring to secure the catheter and avoid kinking the body of the device.

17. The retention clip may be used to secure the catheter or tubing to the patient or drapes. The clip may be initially placed using adhesive tape. To ensure that the catheter is secure during the procedure, clamp or suture the retention clip through the retention clip holes.

Warning: Do not rely on the adhesive alone to secure the clip.

18. When it is determined that extracorporeal bypass is no longer required, clamp the tube section that is not reinforced by wire, carefully remove the catheters from the vessels, and repair their placement sites according to standard surgical practice.

19. After use, this product may be a potential biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable local, regional, and national laws and regulations.

10. The following disclaimer of warranty applies to United States customers only:

Disclaimer of warranty

ALTHOUGH THE CATHETER AND INTRODUCER, HEREAFTER REFERRED TO AS "PRODUCT," HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE, OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT, OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable, or in conflict with applicable law, by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

11. The following disclaimer of warranty applies to customers outside the United States:

Disclaimer of warranty

Important: This disclaimer of warranty does not apply in Australia.

ALTHOUGH THE CATHETER AND INTRODUCER, HEREAFTER REFERRED TO AS "PRODUCT," HAS BEEN CAREFULLY DESIGNED, MANUFACTURED, AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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Saumya Vora Digitally signed by Saumya Vora
Date: 2021.08.26 09:56:08 07'00"

Saumya Vora
Regulatory Affairs Specialist,
Medtronic, Inc.

2019/NOV/13 at 4:28 p.m. Doc number: M871130A001 Rev AA [Ifudoc]
Printspec c - 5.5 x 8.5 inches

Medtronic

Bio-Medicus Life Support™

Catheter and Introducer - Cardiopulmonary Bypass
Katéter a zaváděč - Mímotělní oběh
Kateter og innfører - Kardiopulmonal bypass
Katheter und Einführungshilfe - Herz-Lungen-Maschine
Καθετήρας και εισαγωγικός - Καρδιοπνευμονική παράκαμψη
Catéter e introductor - Bypass cardiopulmonar
Kaletri ja sisäännviejä - sydän-keuhkokoneen käyttö
Cathéter et introducteur - Circulation extracorporelle
Katéter és bevezető - kardiopulmonális bypass
Catetere e introduttore - Bypass cardiopolmonare
Katheter en introducer - Cardiopulmonale bypass
Kateter og innfører - kardiopulmonal bypass
Cewnik i introduktor - krążenie pozaustrojowe
Cateter e introdutor - Circulação extracorpórea
Cateter e introdutor - bypass cardiopulmonar
Катетер и интродьюсер — экстракорпоральное кровообращение
Kateter och introducer - kardiopulmonell bypass
Kateter ve İntrodüser - Kardiopulmoner Bypass


Instructions for Use • Návod k použití • Brugsanvisning • Gebrauchsanweisung • Οδηγός χρήσης •
Instrucciones de uso • Käyttöohjeet • Mode d'emploi • Használati útmutató • Istruzioni per l'uso •
Gebruiksaanwijzing • Bruksanvisning • Instrukcja użytkowania • Instruções de utilização •
Instruções de utilização • инструкция по применению • Bruksanvisning • Kullanım Talimatları

 Caution: Federal law (USA) restricts this device to sale by or on the order of
a physician.

CE
0123

master:20190404, sldocMod:20180927, sldocMod2:20171221, sldocStyle_graph:20190719, sldocStyle:20180308,
templates_global:20180215, master_page_set:20180222, page_specs:20180128, body_ext:20171025, languageLookUp:
20171025, table:20180906, characterJoin:20180903, JoF:20110211, JoF_svg:20110211

Saadat M. Ali 11 Nov 2021

As a registered Notary Public for the State of Minnesota,
I hereby certify that this is a true and accurate copy of
the original document. The original is filed at
Medtronic, Inc.

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Bio-Medicus Life Support™

Catheter and Introducer - Cardiopulmonary Bypass

1. Product description

The Bio-Medicus Life Support catheter is a single-lumen catheter used to drain or reinfuse blood. These devices are sterile, nonpyrogenic, disposable, and intended for single use only. Do not store the product above 25°C (77°F).

2. Indications for use

Bio-Medicus Life Support catheters are peripheral-access catheters used to perfuse vessels or organs in conjunction with extracorporeal cardiopulmonary bypass (CPB) procedures. The catheter introducer is intended to facilitate proper insertion and placement of the appropriately-sized catheter within the vessel. Catheter models with tip lengths of 18 cm (7.09 in), 50 cm (19.7 in), or 55 cm (21.7 in) and without additional side holes may be used as either drainage or reinfusion catheters. This product is intended for use in adult and pediatric patients for up to 6 hours.

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4. Warnings and precautions

- Read all warnings, precautions, and instructions for use carefully before use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.
- Verify the use-by date.
- This device was designed for single use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Do not autoclave.
- Do not use if the sterile package has been damaged or if the seal has been broken in any way.
- Do not use alcohol or alcohol-based fluids for lubrication.
- Thoroughly inspect the catheter before use to verify that the lumen and side holes (if applicable) are patent and that the catheter has not been damaged or kinked.
- Use caution when assembling or manipulating the introducer within the catheter, or moving the suture ring.
- Avoid false lumen catheterization.
- Never advance, position, or place the catheter without the introducer. This action may kink and damage the catheter.
- Use care to protect the side hole areas from kinking as they are passed into the vessel.
- Manipulate the catheter to prevent kinks or any restrictions that may alter flow.
- Manipulation of the heart with the catheter in place may cause flow restrictions.
- Avoid kinking when attaching tubing and when anchoring or untethering the catheter body.
- Use caution when attaching the catheter to other devices in the extracorporeal circuit to prevent damage to the device being connected or the catheter.
- Clamping too close to the connector may damage the connector end.
- Do not clamp wire-reinforced sections. Clamping on reinforced sections may result in permanent wall distortion or lumen collapse.
- Since the catheter wall is opaque due to the presence of the spring, be sure to completely remove the air from the catheter during the priming procedure.
- Avoid overinsertion of the catheter and impingement of the catheter tip against the vessel, ventricular, or atrial walls.
- The use of cold solutions and cardioplegic agents, administered externally or via the coronary circulation, may increase the stiffness of urethane plastic. Increased stiffness may contribute to ventricular perforation under conditions of catheter tip pressure during manipulation of the heart.
- Do not attempt to cross-clamp the catheter with the introducer in place.
- Carefully monitor for signs of occlusion on both the inlet and outlet catheters.
- The user assumes responsibility for any variations in extracorporeal procedures that could compromise the intended use of this device.
- Initiate adequate anticoagulation therapy before patient catheterization.
- Only physicians trained and experienced in using percutaneous catheterization techniques (such as the Seldinger technique) and CPB should use this device.

Instructions for Use English 9

master:20190404, alldocMod:20180927, alldocMod2:20171221, alldocStyle_graph:20190719, alldocStyle:20180308, templates_global:20180215, master_pageset:20180222, pagespecs:20180129, body_extnt:20171025, languageLookUp:20171025, table:20180606, characterJoin:20180903, 3of9:20110211, 3of9_svg:20110211

- Ischemia may occur distal to the catheter placement site.
- Catheter tip placement should be confirmed by using fluoroscopy.
- Use extreme caution to ensure that the return catheter and drainage catheter are not confused while being connected to the perfusion circuit.
- In some cases, suture rings may not be detectable on x-ray due to the variances in anatomical densities and x-ray techniques. Modifying the x-ray angle or distance may aid in detection.
- When utilizing vacuum for augmented drainage, do not exceed -100 mmHg.
- Maximum pressure is 300 mmHg.
- Do not expose to MRI while the catheter is inserted in the patient.

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6. Complications

Possible complications include those normally associated with CPB, anticoagulation, and large bore venous and arterial catheterization.

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1. Select a catheter size that is appropriate for the patient based on criteria including desired flow rate, body surface area, vasculature size, and catheter flow charts. Various sizes are available in order to select the appropriate catheter for each specific procedure.

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3. Using aseptic technique, open the protective, sterile package containing the catheter and introducer.

Note: Use a Bio-Medicus™ insertion kit when performing Seldinger or percutaneous catheterization technique.

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7. Position the adjustable radiopaque suture ring to the desired depth.
8. Pass the catheter and introducer assembly over the guidewire, advancing until the guidewire can be grasped beyond the introducer hub.
9. Insert tip of the catheter introducer into the vessel approximately 3 cm (1.18 in) to 5 cm (1.97 in). Grasp the guidewire and hold stationary while advancing the catheter to the desired depth. Blood leakage from the proximal drainage holes may occur during catheter insertion until these holes are fully advanced into the vessel.
Note: Advancing the catheter and introducer with a rotating motion will make catheter placement easier.
Note: When the bold band closest to the introducer hub is in alignment with the edge of the hemostasis cap, the introducer tip and catheter tip are at the same depth (Refer to Figure 1). For drainage catheter sizes 15 Fr to 29 Fr (5.0 mm to 9.7 mm), the alignment of the catheter and introducer ends is indicated by the tip of the triangle-shaped symbol.
10. **For drainage catheter sizes 21 Fr to 29 Fr (7.0 mm to 9.7 mm) only:** When advancing the venous catheter into the superior vena cava (SVC), verify positioning with tactile feel, indirect visualization method, or echocardiographic imaging.
Note: For drainage models, graduated marks on the hub end of the introducer allow users to approximate the introducer length beyond the tip of the catheter.
11. For proper intracardiac advancement, once the side holes are intravascular, withdraw the introducer to the base of the triangle-shaped symbol representing the tapered tip of the introducer.
Warning: Do not withdraw to the point of the triangle-shaped symbol because of the risk of SVC, innominate vein perforation or vascular, cardiac trauma during intravascular advancement. Never force the catheter. It should advance easily.
12. When desired catheter placement has been achieved, remove the guidewire and retract the introducer clear of the clamp site.
Caution: To maintain hemostasis, do not completely remove the introducer or hemostasis cap until the clamp is applied.
Warning: Ensure that the clamp is properly positioned on the catheter tubing between the tubing connector and wire-reinforced sections. Clamping too close to the connector may damage the connector end. Clamping on wire-reinforced sections may result in permanent wall distortion or lumen collapse.
Warning: When the bold band closest to the vascular end of the introducer is visible at the edge of the hemostasis cap, the catheter can be safely cross-clamped. Do not attempt to cross-clamp the introducer.
Note: Perfusion tubing clamps must be used. Many clamps are too weak to ensure complete clamping.
13. Remove the hemostasis cap, exposing the barbed connector.
Note: Drainage catheter caps are color-coded blue.
Note: The 18 cm (7.09 in), 50 cm (19.7 in), and 55 cm (21.7 in) insertion-length catheters (with side holes at the tip only) may be used as either a return catheter or a drainage catheter. The hemostasis caps for these catheters are clear, not color-coded.
14. Slowly fill the remainder of the catheter tube with sterile priming fluid, displacing any air out through the tube connector end before connecting fluid-filled tubing.
15. Attach the circuit, taking care to connect the drainage catheter to the drainage line (inlet to pump) and the return catheter to the return line (outlet from oxygenator) to ensure an air-free connection. The tubing may be further secured with cable ties at the barbed connector.
Note: The drainage catheter is identified by an arrow in the clamp area next to the Bio-Medicus Life Support logo, indicating flow out of the connector. See Figure 1.
16. Secure the device to the patient following standard medical practice.
Note: In order to avoid accidental displacement during extracorporeal circulation, suture may be wrapped around the suture ring to secure the catheter and avoid kinking the body of the device.
17. The retention clip may be used to secure the catheter or tubing to the patient or drapes. The clip may be initially placed using adhesive tape. To ensure that the catheter is secure during the procedure, clamp or suture the retention clip through the retention clip holes.
Warning: Do not rely on the adhesive alone to secure the clip.

18. When it is determined that extracorporeal bypass is no longer required, clamp the tube section that is not reinforced by wire, carefully remove the catheters from the vessels, and repair their placement sites according to standard surgical practice.
19. After use, this product may be a potential biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable local, regional, and national laws and regulations.

9. The following disclaimer of warranty applies to United States customers only:

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10. The following disclaimer of warranty applies to customers outside the United States:

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©2020 Medtronic. All rights reserved. Medtronic and Medtronic logo are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.