Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA

ATTESTATION STATEMENT


As a responsible member of Medtronic, Inc., I hereby attest that I believe the following Instructions for Use, attached herein, to be true and accurate copies of the original documents.

Alan Gardiner

Alan Gardiner
Principal Regulatory Affairs Specialist

29 Nov 2021

Date

Kathrine Mary Aaberg
November 29, 2021

Medtronic

Alligator Clip Screening Cable
Accessory Kit

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USA Rx only

Medtronic

Manufacturer
Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA

Internet: www.medtronic.com
Tel: 1-763-505-5000
Fax: 1-763-505-1000

Medtronic E.C. Authorized
Representative/Distributed by
Medtronic B.V.
Earl Bakkenstraat 10
5422 PJ Heerlen
The Netherlands
Tel: 31-45-565-8000
Fax: 31-45-565-8668

Europe/Africa/Middle East
Headquarters
Medtronic International Trading Sàrl
Route du Molliau 31
Case Postale 84
CH-1131 Tolochenaz
Switzerland
Internet: www.medtronic.com
Tel: 41-21-802-7000
Fax: 41-21-802-7900

Asia-Pacific
Medtronic International Ltd.
Suite 1602 16/F, Manulife Plaza
The Lee Gardens, 33 Hysan Avenue
Causeway Bay
Hong Kong
Tel: 852-2891-4068
Fax: 852-2591-0313

© Medtronic, Inc. 2010
All Rights Reserved
M928984A007

2011-01

M928984A007 Rev A

Printing instructions. Refer to the
specific Printing Instructions
in doc # A12534.

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Lot number

Număr lot

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

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Lot number

Număr lot

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.

Conformité Européenne
(European Conformity). This
symbol means that the device fully
complies with European Directive
MDD 93/42/EEC.</

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
Rx only

M928983A005 Rev A

Printing Instructions: Refer to the "Addendum" category Table 1 in doc# A00002 for Neuro Core European Printing Instructions.



2010-
January

Filename Date Time
UC200xxxxx EN
4.625" x 8" (117.475 mm x 152.4 mm)

Medtronic Confidential
IPU@sheetmanual.xsl -
IPUGlobalManual@mail.fm
Version: 04-20-2011


M926983A005 Rev A

Printing instructions: Refer to the "Addendum" category Table 1
in doc# A00002 for Neuro Core European Printing Instructions.

2010-
January

Filename Date Time
UC200xxxxx EN
4.625" x 6" (117 475 mm x 152 4 mm)

Medtronic Confidential
IFUSheetManual.xml -
IFUGlobalManualSmall.fm
Version: 04-20-2011

Explanation of symbols on package labeling

Do not reuse

Sterilized using ethylene oxide

Caution, consult accompanying documents

Consult instructions for use

Date of manufacture

Manufacturer

Use by

Temperature limitation

Lot number

Conformité Européenne (European Conformity). This symbol means that the device fully complies with European Directive MDD 93/42/EEC.

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For USA audiences only

Twist-lock Screening Cable Accessory Kit 3550-68 English 1

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Description

The Medtronic® Model 3550-68 Twist-lock Screening Cable Accessory Kit is designed for use with compatible Medtronic Neuromodulation DBS lead kits.

This accessory kit contains one screening cable and one short stylet. The connector insulation is nylon/polycarbonate and the cap insulation is polycarbonate. The short stylet handle material is acetal resin and the wire is tungsten.

The cable is for single use only.

Indications for Use

Refer to the indications sheet packaged with the neurostimulator for indications and related information.

Information for prescribers

Refer to the documentation packaged with the implanted product for important device therapy information.

Warnings

Single use only - This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

Precautions

Storage Temperature - Do not expose the accessories to storage temperatures above 57 °C (135 °F) or below -34 °C (-30 °F).

Instructions for use

Note: Before opening the package, verify the model number and use-by-date.

The twist-lock cable is quadripolar.

⚠ **Warning:** To prevent possible uncomfortable or unexpected stimulation (jolting or shocking sensation), always verify the amplitude setting for the external neurostimulator is 0.0 V before attaching the twist-lock connector to the lead contacts.

⚠ **Caution:** Before connecting the twist-lock screening cable to the stylet handle of the lead, secure the cable to the frame or another stable object. If the cable is not secured, the weight of the twist-lock connector may cause lead movement.

Twist-lock Screening Cable Accessory Kit 3550-68 English 19

M928983A005 Rev A Printing instructions: Refer to the "Addendum" category Table 1 in doc# A00002 for Neuro Core European Printing Instructions. 2010-January

△ **Caution:** Do not pull on the cable. Pulling on the cable may break a wire or dislodge the lead. A broken wire or dislodged lead may result in loss of stimulation and may require surgery to replace the lead.

1. Insert and lock the stylet handle on the lead into the twist-lock connector on the screening cable (Figure 1).

Note: The handle fits into the twist-lock connector in only one way (Figure 1).


Figure 1. Insert the stylet handle into the twist-lock connector.

Note: Place the stylet handle into the groove at a slight angle to secure the handle end (Figure 2).


Figure 2. Place the stylet handle in the groove.

20 English Twist-lock Screening Cable Accessory Kit 3550-68

M928983A005 Rev A Printing Instructions: Refer to the "Addendum" category Table 1 in doc# A00002 for Neuro Core European Printing Instructions.

2010-
January


Figure 3. Press the stylet handle end in the groove.


Figure 4. Turn to lock the twist-lock connector.

2. Ensure the external neurostimulator (ENS) amplitude is set to 0.0V.
3. Push the plug end of the connector cable in the output jack of the ENS. Match the keyed slots of the cable plug and the output jack (Figure 5).


Figure 5. Plug the connector cable into the ENS.

4. Perform test stimulation using the 8870 software applicable to the neurostimulator being implanted. See instructions for test stimulation in the programming guide that accompanies the software.

5. When test stimulation is complete, unlock the twist-lock connector and remove the stylet handle from the lead.
 - a. Hold the test stimulator end of the twist-lock connector stationary in the left hand and turn the lead end of the twist-lock connector counterclockwise with the right hand until the grooves on each side are lined up (Figure 6).


Figure 6. Unlock the twist-lock connector.

- b. Gently pull up on the lead end of the stylet handle until it is free to remove it from the twist-lock connector.
6. After the stylet has been removed and the lead is stabilized, recheck the stimulation effect by reconnecting the twist-lock connector to the lead using the short stylet. Before reconnecting the cable, set the ENS amplitude to 0.0 V.
 - a. Attach the short stylet to the proximal end of the lead (Figure 7).


Figure 7. Attach the short stylet to the lead.

- b. Insert the short stylet completely into the lead (Figure 7a).
 - c. Secure the lead in the stylet handle (Figure 7b).
 - d. Insert the stylet handle of the lead into the twist-lock connector on the screening cable and lock it (refer to Figure 1 through Figure 4).
7. Remove the twist-lock cable from the lead when the recheck is complete.
8. Turn the ENS off, and disconnect the cable from the ENS.

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Medtronic

Alleviating Pain • Restoring Health • Extending Life

Manufacturer and
Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA
Internet: www.medtronic.com
Tel. 1-763-505-6000
Fax 1-763-505-1000

Medtronic E.C. Authorized EC REP
Representative/Distributed by
Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Tel. 31-45-566-8000
Fax 31-45-566-8668

Europe/Africa/Middle East
Headquarters
Medtronic International Trading Sàrl
Route du Molliet 31
Case Postale 84
CH-1131 Tolochenaz
Switzerland
Internet: www.medtronic.com
Tel. 41-21-802-7000
Fax 41-21-802-7900

Asia-Pacific
Medtronic International Ltd.
Suite 1802 16/F, Manulife Plaza
The Lee Gardens, 33 Hysan Avenue
Causeway Bay
Hong Kong
Tel. 852-2591-4068
Fax 852-2591-0313

* M 9 2 8 9 8 3 A 0 0 5 *

© Medtronic, Inc. 2010
All Rights Reserved
M928983A005

M928983A005 Rev A

2010-
January

Filename Date Time
UC200xxxxx EN
4.625 x 6 inches (117 mm x 152 mm)

Medtronic Confidential
Indications.xsl - IndicationsTemplate.fm
Template Version: 10-06-2010

Medtronic


Medtronic® DBS™ Therapy for Movement Disorders
Thérapie DBS™ de Medtronic® pour le traitement des troubles de la motricité
Medtronic® DBS™-Therapie für Bewegungsstörungen
Terapia DBS™ de Medtronic® para trastornos de movimiento
Medtronic® DBS™-therapie voor bewegingsstoornissen
Terapia DBS™ della Medtronic® per il trattamento dei disturbi del movimento
Medtronic® DBS™-behandling for bevægelsesforstyrrelser
Medtronic® DBS™ -hoito liikehäiriöiden hoitoon
Medtronic® DBS™-behandling för rörelserubbningar
Θεραπεία DBS™ της Medtronic® για κινητικές διαταραχές
Medtronic® DBS™-behandling for bevægelsesrelatede lidelser
Terapia DBS™ da Medtronic® para a epilepsia
Medtronic® DBS™ terapia mozgászavarok kezelésére
Terapia Medtronic® DBS™ terapia do leczenia zaburzeń ruchowych
Systém DBS™ společnosti Medtronic® pro léčbu poruch pohybového aparátu
Терапия двигательных расстройств с помощью системы
Medtronic® DBS™
DBS™ na liečbu porúch pohybového aparátu od spoločnosti Medtronic®
Hareket Bozukluklarına yönelik Medtronic® DBS™ Tedavisi
Terapia DBS™ pentru tulburări de mișcare Medtronic®
علاج Medtronic® DBS™ لأمراض اضطرابات الحركة
טיפול Medtronic® DBS™ בהפרעות תנועה

Indications sheet • Feuille d'indications • Indikationsblatt •
Hoja de indicaciones • Indicatieboekje • Foglietto illustrativo •
Vedlegg om indikationer • Käyttöaihelehtinen • Lista över indikationer •
Φύλλο ενδείξεων • Indikationsside • Folha de indicações •
Javallatokat ismertető lap • Informacje o wskazaniach •
Informace o indikacích • Перечень показаний • Zoznam indikácií •
ورقة دواعي الاستعمال • גיליון התוויות • Listă de indicații

M930110A002 Rev A 2010-10



Printing instructions: Refer to the "Addendum" category
Table 1 in doc# A00002 for Neuro Core European
Printing Instructions.

Filename Date Time
UC200xxxxxx EN
4.625 x 6 inches (117 mm x 152 mm)


Medtronic Confidential
Indications.xsl - IndicationsTemplate.fm
Template Version: 10-06-2010

M930110A002 Rev A 2010-10

Printing instructions: Refer to the "Addendum" category
Table 1 in doc# A00002 for Neuro Core European
Printing Instructions.

Filename Date Time
UC200xxxxxx EN
4.625 x 6 inches (117 mm x 152 mm)

Medtronic Confidential
Indications.xsl - IndicationsTemplate.fm
Template Version: 10-06-2010

Medtronic® is a registered trademark of Medtronic, Inc.
DBS™ is a trademark of Medtronic, Inc.

Indications English 1

M930110A002 Rev A 2010-10

Filename Date Time
UC200xxxxx EN
4.625 x 6 inches (117 mm x 152 mm)

Medtronic Confidential
Indications.xsl - IndicationsTemplate.fm
Template Version: 10-06-2010

Indications

Medtronic DBS Therapy for Movement Disorders is Indicated for patients with disabling tremor or symptoms of Parkinson's disease. Studies have shown that deep brain stimulation with Medtronic DBS Therapy components is effective in controlling essential tremor and symptoms of Parkinson's disease that are not adequately controlled with medications. Additionally, deep brain stimulation is effective in controlling dyskinesias and fluctuations associated with medical therapy for Parkinson's disease.


Medtronic DBS Therapy for Movement Disorders is also indicated as an aid in the management of chronic, intractable (drug refractory) primary dystonia, including generalized and segmental dystonia, hemidystonia, and cervical dystonia (torticollis) for individuals 7 years of age and older.

Clinicians, refer to the appropriate information for prescribers booklet for contraindications, warnings, precautions, adverse events summary, patient selection, and component disposal.

Patients, refer to your Medtronic DBS patient therapy guide for contraindications, warnings, precautions, and other risk information.

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








Medtronic
Alleviating Pain • Restoring Health • Extending Life

Manufacturer
Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA
Internet: www.medtronic.com
Tel. 1-763-505-5000
Fax 1-763-505-1000

Medtronic E.C. Authorized **EC** **REP**
Representative/Distributed by
Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Tel. 31-45-566-8000
Fax 31-45-566-8068

Europe/Africa/Middle East
Headquarters
Medtronic International Trading Sàrl
Route du Molkau 31
Case Postale 84
CH-1131 Tolochenaz
Switzerland
Internet: www.medtronic.com
Tel. 41-21-802-7000
Fax 41-21-802-7900

Asia-Pacific
Medtronic International Ltd.
Suite 1602 16/F, Manulife Plaza
The Lee Gardens, 33 Hysan Avenue
Causeway Bay
Hong Kong
Tel. 852-2891-4068
Fax 852-2591-0313

© Medtronic, Inc. 2010
All Rights Reserved
M930110A002

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! USA Rx only

CE0123

Medtronic® is a trademark of Medtronic, Inc., registered in the U.S. and other countries.

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USA FCC Information

The following is communications regulation information on the Model 37021 and 37022 external neurostimulators.

FCC ID: LF537021

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

IMPORTANT: Changes or modifications to this product not authorized by Medtronic, Inc., could void the FCC Certification and negate your authority to operate this product.

Contents

Label symbols 5

1 Introduction 7

Purpose of the device 8

Precautions 9

How to use this guide 10

Test stimulation overview 10

Test stimulation goals 11

Patient guides 12

2 Using your test stimulation system 13

Knowing your system components 14

Using your external neurostimulator and accessories 16

Using the external neurostimulator holder 16

Checking the external neurostimulator batteries 18

Replacing the external neurostimulator batteries 19

Using the Therapy-stop key 22

Reconnecting the cable to the external neurostimulator 22

Specifications 24

Safety and technical checks 26

Battery and device disposal 27

Declaration of conformity 28

3 Your daily activities 29

What you should know about your activities 30

Device care 30

Physical activities 31

Incision care 32

Sleeping 33

What changes in sensation to expect during test
stimulation 33

Index 35

Contents

4 English 37022 2013-11

Label symbols

Explanation of symbols on products and packaging.
Refer to the appropriate product to see symbols that apply.

CE0123 Conformité Européenne (European Conformity).
This symbol means that the device fully complies with MDD 93/42/EEC (NB 0123) and R&TTE Directive 1999/5/EC.

Consult instructions for use

IEC 60601-1/EN60601-1, Type BF Equipment

Non-ionizing electromagnetic radiation

System meets the applicable Canadian (CAN/CSA-C22.2 No. 60601-1) electrical safety standard requirements.

For USA audiences only

Temperature limitation

Ingress protection rating IP22, per 60601-1-11

Keep dry

Serial number

PIN number

Label symbols

Chinese Standard (SJ/T11364-2006) Logo:
Electronic Information Products Pollution Control
Symbol. (The date in this logo means the
environmental protection use period of the product.)

Do not dispose of this product in the unsorted
municipal waste stream. Dispose of this product
according to local regulations. See [http://
recycling.medtronic.com](http://recycling.medtronic.com) for instructions on proper
disposal of this product.

Manufacturer

Authorized representative in the European
community

Magnetic Resonance (MR) Unsafe

Label symbols

6 English 37022 2013-11

37022_Ch01.fm 1/31/14 9:52 am
Size: 4.625" x 6.0" (117 mm x 152 mm)
JC 00601-1-11 Update

Medtronic Confidential
NeuroPain®R00

1 Introduction

M520468A023 Rev A

Purpose of the device

The Model 37022 external neurostimulator provides information for evaluating the success of Medtronic neurostimulation systems for pain therapy prior to implant.

The Medtronic neurostimulation system is intended to provide electrical stimulation to relieve chronic, intractable pain. The Medtronic neurostimulation system for pain therapy indications, contraindications, risks and benefits, warnings, precautions, and individualization of treatment information are located in the *Patient Programmer: Pain Therapy User Manual*.

Precautions

Equipment modification – Do not modify this equipment. Modification of this equipment can result in damage to the device causing the device to malfunction or become unusable.

Equipment use – When operating an external neurostimulator, use special care near flammable or explosive atmospheres. An interaction between the flammable or explosive atmospheres and the battery in the device could occur. The consequences of using a battery-powered device near flammable or explosive atmospheres are unknown.

Interference from EMI – Do not use the device in the proximity of equipment that generates electromagnetic interference (EMI). EMI may cause a disruption in device function. Examples of common medical sources of EMI are magnetic resonance imaging (MRI) and lithotripsy. Powerful computer monitors, cell phones, x-ray equipment, and other monitoring equipment may also generate EMI.

How to use this guide

For over 30 years, stimulation has helped thousands of patients manage their pain and improve their quality of life. Test stimulation will determine if stimulation is helpful for you. Stimulation will not cure your pain. It may, however, reduce your pain to a tolerable level. As a result, you may be able to resume your daily activities.

This guide includes the following information:

- Chapter 1, "Introduction," provides an overview and goals of test stimulation.
- Chapter 2, "Using your test stimulation system," describes how to use your external neurostimulator.
- Chapter 3, "Your daily activities," gives recommendations regarding your daily activities during test stimulation.

Ask your doctor to explain anything that is unclear.

Test stimulation overview

Your pain is unique to you. Because pain relief can only be evaluated by you, your doctor has prescribed a test stimulation period. Test stimulation is designed to help you and your doctor evaluate the benefits of stimulation for managing your pain.

Some doctors may prescribe three to five days for test stimulation; others may prescribe more days. Your doctor will determine the number of days that is best for you.

Test stimulation allows you to experience stimulation during your daily activities. Your doctor may ask you to keep a log of how you feel at different times during the day. The log helps your doctor determine if the stimulation is effective for managing your pain.

If you experience adequate pain relief during test stimulation, and you are comfortable using the system, your doctor may decide to implant a neurostimulation system.

Test stimulation goals

The goals of a test stimulation are to determine:

- If stimulation covers your body's painful areas
- If you are comfortable with the sensation of stimulation
- If you can use the patient programmer to control your stimulation

Patient guides


In addition to the information contained in this guide, your clinician should have provided you with the following patient materials:

- ***Patient Programmer: Pain Therapy User Manual.***
This manual provides an overview of how stimulation works, a detailed description of how to use your patient programmer, and troubleshooting steps.
- ***Patient Programmer: Quick Reference Guide.***
This guide provides basic, step-by-step instructions on how to use your patient programmer.

Review these materials and keep them in a convenient place during test stimulation.

37022_CN02.fm 10/1/14 9:52 am
Size 4 625" x 6 0" (157 mm x 152 mm)
UC 60601-1-11 Update

Medtronic Confidential
NeuroPainIR00

2 Using your test stimulation system

M928468A023 Rev A

Knowing your system components

For test stimulation, your clinician will provide you with the following components:

- External neurostimulator
 - Two AA alkaline batteries (LR6)^a
- Connector cable
- Patient programmer
 - Two AAA alkaline batteries (LR03)^a
- External neurostimulator holder with belt clip

Using your test stimulation system 2

^a You may need to replace batteries during test stimulation; have extra batteries available.


Figure 2.1 Test stimulation system components.

Using your external neurostimulator and accessories

This section provides you with the instructions for:

- Using the external neurostimulator holder
- Checking the external neurostimulator batteries
- Replacing the external neurostimulator batteries
- Using the Therapy-stop key
- Reconnecting the cable to the external neurostimulator

Note: Before placing the external neurostimulator into operation, ensure the external neurostimulator has had time to equalize to the current temperature and environment.

Using your test stimulation system 2

Using the external neurostimulator holder

Use the external neurostimulator holder to carry the external neurostimulator on your belt or garment.

1. Slide the external neurostimulator holder over your belt or garment (Figure 2.2).

Note: A belt with device pouches is available as an option. Refer to the instructions provided with the belt for additional information.


Figure 2.2 Slide the external neurostimulator holder over your belt or garment.

2. With the triangles on the cable plug and the external neurostimulator pointing toward each other (Figure 2.3), gently insert the plug fully into the output jack.


Figure 2.3 Insert the plug into the output jack.

3. Place the external neurostimulator into the holder with the red Therapy-stop key facing out (Figure 2.4).


Figure 2.4 Place the external neurostimulator into the holder.

Checking the external neurostimulator batteries

Use your patient programmer to check your external neurostimulator battery level every day. For instructions on checking the external neurostimulator battery level, refer to the manual packaged with the patient programmer.

If the external neurostimulator batteries need immediate replacement, the patient programmer will display an information or warning screen. Information and warning screens are described in the Troubleshooting section of the manual packaged with your patient programmer.

If test stimulation lasts seven days or less, you should not need to replace your batteries. Battery usage

depends on your specific therapy settings and length of test stimulation.

Replacing the external neurostimulator batteries

Cautions:

- When replacing batteries during test stimulation, insert the new batteries within 15 minutes of removing the old batteries. If the batteries are not replaced within 15 minutes, the neurostimulator may reset. When a neurostimulator resets, stimulation history is no longer available, and the stimulation settings may not reflect recent programming changes.
 - Before installing batteries, always decrease all amplitudes to the lowest setting, turn off the external neurostimulator and then disconnect the cable from the external neurostimulator to prevent possible uncomfortable or unexpected stimulation (jolting or shocking sensation) when stimulation is turned on.
1. Use your patient programmer to decrease all amplitudes and turn off the external neurostimulator.
 2. Remove the external neurostimulator from the holder.
 3. Gently pull the cable plug from the output jack.

Using your test stimulation system 2

⚠ Caution: Do not pull on the cable. Pulling on the cable may break a wire or dislodge the lead. A broken wire or dislodged lead may result in loss of stimulation and may require surgery to replace the lead.

4. Press down lightly on the battery compartment cover (Figure 2.5a).


Figure 2.5 Open the battery compartment cover.

5. Push the cover in the direction of the arrow, then swing the cover open (Figure 2.5b, c).
6. Remove the depleted batteries and insert two new AA alkaline batteries into the battery compartment (Figure 2.5c).

Note: Match the battery polarity to the positive (+) and negative (–) markings on the inside of the battery compartment cover (Figure 2.5c).

7. Push down the battery compartment cover until it clicks into position. Then slide the cover back in the opposite direction of the arrow until it closes.

Using your test stimulation system 2

8. Reconnect the cable. (For detailed instructions, see "Reconnecting the cable to the external neurostimulator" on page 22.)
9. Place the external neurostimulator into the holder.
Note: After the batteries are installed and the battery compartment cover is closed, stimulation may not be available for up to 60 seconds.
10. Use the patient programmer to turn on your external neurostimulator.
11. Dispose of depleted batteries according to local requirements.

Using the Therapy-stop key

Typically, you will use the patient programmer to turn stimulation off. However, if you need to turn stimulation off immediately, press the red Therapy-stop key (Figure 2.6). You must use the patient programmer to turn the external neurostimulator on again.


Figure 2.6 Use the Therapy-stop key to immediately turn stimulation off.

Using your test stimulation system 2

Reconnecting the cable to the external neurostimulator

If the cable plug is disconnected from the external neurostimulator, stimulation will stop.

1. With the triangles on the cable plug and the external neurostimulator pointing toward each other, gently insert the plug fully into the output jack (Figure 2.7).

Caution: Do not pull on the cable. Pulling on the cable may break a wire or dislodge the lead. A broken wire or dislodged lead may result in loss of stimulation and may require surgery to replace the lead.


Figure 2.7 Insert the plug into the output jack.

2. Use your patient programmer to turn stimulation on.

Specifications

Table 2.1 External neurostimulator specifications

Item	Specification
Power source	2 AA alkaline batteries (nonrechargeable, LR6)
Battery life	144 hours minimum (21 days average) for alkaline batteries
System operating temperature ^a	10 °C to 43 °C (50 °F to 109 °F)
Temperature limitation ^{b,c}	-34 °C to 57 °C (-30 °F to 135 °F)
Ingress protection	IP22 rating for solid objects greater than or equal to 12.5 mm, and for vertically dripping water when the device is tilted 15 degrees, per 60601-1-11.
Size (approximate)	8.9 cm x 5.4 cm x 2.8 cm (3.5 in x 2.1 in x 1.1 in)
Weight, including batteries (approximate)	100 g (3.5 oz)
Mode of operation	Continuous
Case material	Polycarbonate/ABS blend plastic resin
External neurostimulator automatically shuts off when ^d :	
• Device detects moisture • Battery door opens • Connector cable is not connected • Connector lid opens	

Using your test stimulation system 2

- ^a Product remains safe when operated at temperatures from 5 °C (41 °F), to less than or equal to 43 °C (109 °F), in relative humidity from 15% to 93%, non-condensing, and in atmospheric pressures from 700 hPa to 1060 hPa.
- ^b Batteries should be removed from the device for storage or transport.
- ^c Product remains safe when stored or transported in temperatures from -34 °C (-30 °F) without relative humidity control, to less than or equal to 70 °C (158 °F) at a relative humidity of less than or equal to 93%, non-condensing.
- ^d Use the patient programmer to turn the external neurostimulator on when the condition is resolved.

Safety and technical checks

Periodic safety and technical checks or periodic maintenance of the external neurostimulator are not required.

The external neurostimulator contains no user-serviceable parts. If repair or service is needed, contact your clinician or a Medtronic sales office. Refer to the list of Medtronic contacts at the end of this manual.

Using your test stimulation system 2

Battery and device disposal

At the end of the test stimulation, return the device to your doctor for disposal. Dispose of depleted batteries according to local requirements. Do not incinerate.

Declaration of conformity

Medtronic declares that this product is in conformity with the essential requirements of Directive 1999/5/EC on Radio and Telecommunications Terminal Equipment and Directive 93/42/EEC on Medical Devices.

For additional information, contact the appropriate Medtronic representative listed on the inside back cover of this manual


Using your test stimulation system 2

28 English 37022 2013-11

M928468A023 Rev A

37022_Ch03.fm 1/31/14 9:52 am
Size 4.625" x 6.0" (117 mm x 152 mm)
UC 60501-1.11 Update

Medtronic Confidential
NeuroPatntR00

3 Your daily activities

M928468A023 Rev A

What you should know about your activities

To reduce the risk of moving your implanted lead(s) after surgery, always follow your clinician's instructions during test stimulation.

Read and understand this patient guide and the other patient materials about your system.

Keep a record of your stimulation parameter settings, daily activities, and pain relief.

Device care

- Keep fresh batteries available.
- Use the patient programmer daily to check the external neurostimulator batteries.
- Replace low or depleted batteries.
- Do not pull on the cable or lead.
- Keep the device out of the reach of children and pets. Keep batteries away from children. If swallowed, contact a doctor at once.

Your daily activities 3

- Handle the device and system components with care. Do not drop, strike, or step on the device or system components.
- Do not dismantle or tamper with the device.
- Clean the outside of the device with a damp cloth when necessary. Mild household cleaners will not damage the device or labels.
- The device and system components are not waterproof. Do not allow moisture to get inside the device or system components.

Physical activities

To reduce the risk of moving your implanted lead(s), always follow your clinician's instructions and these guidelines:

- Limit your physical activities to low or moderate levels.
- Avoid reaching over your head.
- Avoid turning side to side.
- Avoid climbing too many stairs.
- Avoid sitting in a chair for a long period of time.
- Avoid sudden, rapid movements such as twisting, stretching, or bending.

Note: You can move your body in bed without twisting by moving your shoulders and hips at the

Your daily activities 3

same time with a "log-rolling" movement. Ask your clinician to teach this movement to you while you are in the hospital.

- Avoid bending forward, backward, or from side to side.
- Avoid lifting more than 2.3 kilograms (5 pounds).
- Talk to your doctor about when it is appropriate for you to resume sexual activity.

Incision care

- Keep your incision site(s) dry by taking sponge baths instead of showers to help avoid infections.
- Follow your clinician's instructions about cleaning the incision site(s) or changing dressings.
- Follow your clinician's advice about managing your incision pain.
- Call your clinician if you notice any of these signs of an infection:
 - Increased pain in your incision
 - Fever
 - Chills
 - Pus
 - Redness

Your daily activities 3

Sleeping

- Do not sleep on your stomach.
- Sleep on your back or on the side opposite the device.
- Use a firm mattress that supports your legs and back equally.

What changes in sensation to expect during test stimulation

The sensation of stimulation is different from person to person. Most people report a steady, tingling sensation in the area of their pain.

Generally, people experience a constant sensation of stimulation. However you may feel some normal stimulation changes when you suddenly move or change position.

For example, for spinal cord stimulation patients, leaning back may cause the lead to move closer to your spinal cord; this can increase your sensation of stimulation. Other movements may cause the lead to move further away from your spinal cord and decrease your sensation of stimulation.

Typically, you will use the patient programmer to turn stimulation off. However, if you feel a sudden and unpleasant change in stimulation, press the red

**Therapy-stop key on your external neurostimulator to
turn off the stimulation immediately.**

**If you use the red Therapy-stop key, you must use the
patient programmer to turn your stimulation on again.**

Your daily activities 3

34 English 37022 2013-11

Index

Activities 31
Batteries, external neurostimulator
 checking 18
 precaution 9
 replacing 19
Battery level 18
Cable, reconnecting to external neurostimulator 22
Caring for your device 30
Caring for your incision 32
Cleaning device 31
Components 14
External neurostimulator
 batteries
 checking 18
 precaution 9
 replacing 19
 cable, reconnecting 22
 description 8
 holder 16
Holder 16
Incision care 32
Patient programmer
 checking external neurostimulator battery 18
 typical use 22
Precautions 9
Repair 26
Sensation of stimulation 33
Sleeping 33
Specifications 24
System components 14
Test stimulation
 goals 11
 overview 10
Therapy-stop key 22
Turning stimulation off immediately 22

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Belgium: Medtronic Belgium S.A.
Tel. 02-456-0900
Fax 02-460-2667

Canada: Medtronic of Canada Ltd.
Tel. +1-905-460-3800
Fax (1905)-826-6620

Czech Republic: Medtronic Czechia s.r.o.
Tel. 2-965-795-80
Fax 2-965-795-89

Denmark: Medtronic Danmark A/S
Tel. 45-32-48-18-00
Fax 45-32-48-18-01

Finland: Medtronic Finland Oy/Ltd
Tel. (09)-755-2500
Fax (09)-755-25018

France: Medtronic France S.A.S.
Tel. 01-5538-1700
Fax 01-5538-1800

Germany: Medtronic GmbH
Tel. (02159)-81490
Fax (02159)-8149100

Greece: Medtronic Hellas S.A.
Tel. 210-67-79-099
Fax 210-67-79-399

Hungary: Medtronic Hungaria Kft.
Tel. 1-889-06-00
Fax 1-889-06-99

Ireland: Medtronic Ireland Ltd.
Tel. (01)-890-6522
Fax (01)-890-7220

Italy: Medtronic Italia SpA.
Tel. 02-241371
Fax 02-241381

Italy: Medtronic Italia SpA
Tel. 06-328141
Fax 06-3215812

Japan: Medtronic Japan
Tel. 3-6430-2016
Fax 3-6430-7110

Latin America: Medtronic, Inc.
Tel. (1305)-500-9328
Fax (1786)-709-4244

The Netherlands: Medtronic B.V.
Tel. (045)-566-8000
Fax (045)-566-8668

Norway: Medtronic Norge AS
Tel. 67-10-32-00
Fax 67-10-32-10

Poland: Medtronic Poland Sp. z o.o.
Tel. (022)-465-69-00
Fax (022)-465-69-17

Portugal: Medtronic Portugal, Lda.
Tel. 21-724-5100
Fax 21-724-5199

Russia: Medtronic Russia
Tel. (8495) 580-7377
Fax (8495) 580-7378

Slovakia: Medtronic Slovakia, o.z.
Tel. 0268 206 911
Fax 0268 206 999

Spain: Medtronic Ibérica, S.A.
Tel. 91-625-0400
Fax 91-650-7410

Sweden: Medtronic AB
Tel. 08-568-585-00
Fax 08-568-585-01

Switzerland: Medtronic (Schweiz) AG
Tel. 031-868-0100
Fax 031-868-0199

Turkey: Medtronic Turkey
Tel. +90 216 636 1000
Fax +90 216 636 1008

U.K.: Medtronic U.K. Ltd.
Tel. 01923-212213
Fax 01923-241004

USA: Medtronic, Inc.
Tel. (1763)-505-5000
Fax (1763)-505-1000
Toll-free (1800)-328-0810

M928468A023 Rev A

37022_bocna 1/31/14 9:54 am
Size 4.525" x 6.0" (117 mm x 152 mm)
UC 60601-1-11 Update

Medtronic Confidential
NeuroPatentR06

Manufacturer
Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432-5804
USA

www.medtronic.com
Tel. +1-763-605-5000
Fax +1-763-605-1000

Authorized Representative 
in the European Community

Medtronic B.V.
Earl Bakkenstraat 10
5422 PJ Heerlen
The Netherlands
Tel. +31-45-596-6000
Fax +31-45-596-6666

Europe/Africa/Middle East
Headquarters
Medtronic International Trading SARL
Route du Mollet 31
Case Postale 84
CH-1131 Tolochenaz
Switzerland

www.medtronic.eu
Tel. +41-21-902-7000
Fax +41-21-902-7900

Asia-Pacific
Medtronic International Ltd
Suite 1105-11 11/F, Tower 1, The Gateway
25 Canton Road, Tsimshatsui
Kowloon
Hong Kong
Tel. +852-2919-1300
Fax +852-2991-6830

Contacts for specific countries are listed inside this
cover

* K 9 2 8 4 6 8 A 0 2 3 *

© Medtronic, Inc. 2010
All Rights Reserved
M020468A023 Rev A

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
CE0123

M928467A007 A



2013-11
Printing instructions: Refer to the "Implant Manual"
category Table 1 in doc# A00002 for Neuro Core
European Printing Instructions.

Filename Date Time
UC200xxxxx EN
4.625 x 6 inches (117 mm x 152 mm)


Medtronic Confidential
ImplantManual.xsl - IPGTemplate.fm
Template version 6.1: 01-09-2014

M928487A007 A

2013-11
Printing instructions: Refer to the "Implant Manual"
category Table 1 in doc# A00002 for Neuro Core
European Printing Instructions.

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CE0123 Conformité Européenne (European Conformity). This symbol means that the device fully complies with MDD 93/42/EEC (NB 0123) and R&TTE Directive 1999/5/EC

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

Chinese Standard (SJ/T11364-2006) Logo: Electronic Information Products Pollution Control Symbol. (The date in this logo means the environmental protection use period of the product.)

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
Medtronic® is a trademark of Medtronic, Inc., registered in the US and other countries.

[USA] FCC Information

The following is communications regulation information on the Model 37021 and 37022 External Neurostimulators.

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


M928467A007 A


Table of contents


Purpose of the device	9
Description	9
Package contents	9
Accessories	9
Device specifications	9
Declaration of Conformity	13
Instructions for use	14
Replacing the external neurostimulator batteries	14
Attaching the external neurostimulator and connector cable to the clinician programmer	15
Using the external neurostimulator during test stimulation	15
Device care and storage	16
Safety and technical checks	16

Refer to the indications sheet for indications and related information.

Refer to the appropriate information for prescribers booklet for contraindications, warnings, precautions, adverse events summary, individualization of treatment, patient selection, use in specific populations, resterilization, and component disposal.

Refer to System Eligibility, Battery Longevity, Specifications reference manual for neurostimulator selection, battery longevity calculations and specific neurostimulator specifications.


USA Refer to the clinical summary booklet for information on the clinical study results of the neurostimulation system and individualization of treatment.


Filename Date Time
UC200xxxxxx EN
4.625 x 6 inches (117 mm x 152 mm)

Medtronic Confidential
ImplantManual.xsl - IPGTemplate.fm
Template version 6.1: 01-09-2014


6 English 2013-11

2013-11

M928467A007 A

Purpose of the device

The Medtronic Model 37022 External Neurostimulator (ENS) is used to evaluate a Medtronic Neurostimulation System during lead placement or test stimulation. The external neurostimulator and the implantable neurostimulator produce comparable symptom suppression when set to the same parameter settings.

Description

The Medtronic Model 37022 External Neurostimulator is part of a neurostimulation system.

Package contents

- External neurostimulator
- External neurostimulator holder
- AA alkaline batteries (2)
- Product literature
- [USA] Warranty card

Accessories

- Connector cable (packaged separately)

Device specifications

The external neurostimulator (Figure 1) is a programmable device that delivers stimulation through one or more leads. The stimulation settings are stored in programs. A program is a specific combination of pulse width, rate, and amplitude settings acting on a specific electrode combination. Up to four programs can be combined into a group. When using more than one program, the pulses are delivered sequentially—first a pulse from one program, then a pulse from the next program.

Pulse width, amplitude, and electrode polarity for each program within a group can have different values. Rate, rate limits, SoftStart/Stop value, and cycling for each program within a group have the same values.


Figure 1. Model 37022 external neurostimulator.

Table 1. Operating values for the Model 37022 external neurostimulator^a

Programmable parameter	Operating range and resolution ^b
Number of defined groups	1 to 8 ^c
Number of programs per group	1 to 4 ^c
Electrode configuration	2 to 16 electrodes as anode, cathode, or off
Amplitude	0 to 10.5 V with 0.05-V or 0.1-V resolution 0 to 25.5 mA with 0.1-mA resolution
Amplitude – upper patient limit	Programmed value to 10.5 V (same resolution as amplitude) Programmed value to 25.5 mA (same resolution as amplitude)
Amplitude – lower patient limit	0 V to the programmed value (same resolution as amplitude) 0 mA to the programmed value (same resolution as amplitude)
Pulse width	60 to 1000 μ s (10- μ s resolution)
Pulse width – upper patient limit	Programmed value to 1000 μ s (10- μ s resolution)
Pulse width – lower patient limit	60 μ s to the programmed value (10 μ s resolution)
Rate	2 to 1000 Hz ^d
Rate – upper patient limit	Programmed value to 1000 Hz (same resolution as rate)
Rate – lower patient limit	2 to the programmed value
SoftStart/Stop	Off, on: 1, 2, 4, or 8 second ramp duration

Table 1. Operating values for the Model 37022 external neurostimulator^a (continued)

Programmable parameter	Operating range and resolution ^b
Cycling	Off, on: 0.1 s to 24 hrs
Scheduled therapy	Off, on: 1 to 8 events per 24-h period

- ^a Selection of parameters and values is limited by the clinician programmer with regards to appropriate applications, therapies, and devices
- ^b Interlocks and out-of-regulation detection will prevent the use of some parameter combinations
- ^c No more than 16 programs may be defined within the 8 groups.
- ^d Rate and rate upper limit for the ENS is limited to 500 Hz when two programs are active in a group, 330 Hz when three programs are active in a group and to 250 Hz when four programs are active in a group.

Table 2. Physical characteristics of the Model 37022 external neurostimulator^a

Description	Value
Length	89 mm (3.5 in)
Width	54 mm (2.1 in)
Thickness	28 mm (1.1 in)
Weight (with batteries)	100 g (3.5 oz)
Output jack	21 pin
Identification code	NJS

- ^a All measurements are approximate

Table 3. External neurostimulator electrical and operating characteristics

Description	Value
Power source	AA alkaline batteries (2)
Battery life	144 hours minimum (21 days average) for alkaline batteries ^a
External neurostimulator service life	Up to 5 years
Operating type	Continuous
Degree of protection against electrical shock	Type BF
Ingress protection	IP22 rating for solid objects greater than or equal to 12.5 mm, and for vertically dripping water when the device is tilted 15 degrees, per 60601-1-11.
Case material	Polycarbonate/ABS blend plastic resin
Automatic shut off ^b	- Moisture detected - Battery door open - Connector cable disconnected - Connector lid open

- ^a Battery life is based on 24-hour per day stimulation using 1 group containing 3 active programs. For all programs: impedance = 500 Ω , Amp = 4.0 V, PW = 330 μ s, Rate = 70 Hz, Cycling = off.
- ^b Use the clinician programmer or patient programmer to turn the external neurostimulator on once the condition is resolved.

Table 4. External neurostimulator storage and operating conditions

Parameter	Temperature limitation ^{a, b}	System Operating ^c
Minimum temperature	-34 °C (-30 °F)	10 °C (50 °F)
Maximum temperature	57 °C (135 °F)	43 °C (109 °F)

- ^a Batteries should be removed from the device for storage or transport.
- ^b Product remains safe when stored or transported in temperatures from -34 °C (-30 °F) without relative humidity control, to less than or equal to 70 °C (158 °F) at a relative humidity of less than or equal to 93%, non-condensing.
- ^c Product remains safe when operated at temperatures from 5 °C (41 °F), to less than or equal to 43 °C (109 °F), in relative humidity from 15% to 93%, non-condensing, and in atmospheric pressures from 700 hPa to 1060 hPa.

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For additional information, contact the appropriate Medtronic representative listed on the inside back cover of this manual.

2013-11 English 13

M928467A007 A

2013-11

Instructions for use

The external neurostimulator is used to evaluate lead placement and stimulation settings.

- △ **Caution:** The device is not certified for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide. The consequences of using the device near flammable atmospheres are unknown.
- △ **Caution:** Do not modify this equipment. Modification of this equipment can result in damage to the device, causing the device to malfunction or become unusable.
- △ **Caution:** Do not use the device in the proximity of equipment that generates electromagnetic interference (EMI). EMI may cause a disruption in device function. Examples of common medical sources of EMI are magnetic resonance imaging (MRI) and lithotripsy. Powerful computer monitors, cell phones, x-ray equipment, and other monitoring equipment may also generate EMI.

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Replace the external neurostimulator batteries before each test stimulation and when the batteries are low or depleted. The battery level is shown on the clinician programmer screen.

- △ **Cautions:**
 - When replacing batteries during test stimulation, insert the new batteries within 15 minutes of removing the old batteries. If the batteries are not replaced within 15 minutes, the neurostimulator may reset. When a neurostimulator resets, stimulation history is no longer available, and the stimulation settings may not reflect recent programming changes.
 - Do not leave depleted batteries in the external neurostimulator. The batteries may corrode and cause damage to the electronic components.
 - If the device will not be used for several weeks, remove the batteries from the device. A battery left in the device may corrode, causing damage to the electronic components.
1. If the external neurostimulator is on, use the clinician programmer or patient programmer to turn the neurostimulator off.
 2. If the connector cable is attached to the external neurostimulator, disconnect the cable.
 3. If the neurostimulator is in the external neurostimulator holder, remove the neurostimulator from the holder.
 4. Press down lightly on the battery compartment cover, push the cover in the direction of the arrow, then swing the cover open.
 5. Insert two new AA alkaline batteries. Correct battery polarity is indicated on the battery compartment cover.

6. Replace the battery compartment cover, then slide the cover in the opposite direction of the arrow until it snaps into place.
7. Reconnect the cable. For detailed instructions, see page 15.

Notes:

- After the batteries are installed and the battery compartment cover is closed, the neurostimulator may take up to 60 seconds for device initiation. Stimulation is not available until device initiation is complete.
- Dispose of depleted batteries and worn out devices according to local requirements.

Attaching the external neurostimulator and connector cable to the clinician programmer

The external neurostimulator is usually attached to the clinician programmer during lead placement.

1. Locate the programming head on the back of the clinician programmer.
2. Place the external neurostimulator so it fits into the slot on the programming head.
3. Turn the neurostimulator so the output jack is positioned as displayed in Figure 2.

Note: Facing the output jack away from the clinician programmer properly aligns the two devices for communication and prevents the cable from resting on the programmer.

Figure 2. External neurostimulator attached to the clinician programmer.

Using the external neurostimulator during test stimulation

When programming during test stimulation, hold the clinician programmer programming head over the external neurostimulator. The neurostimulator does not need to be attached to the programmer.

Connecting the connector cable to the external neurostimulator

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- Matching the keyed slots of the cable plug and the output jack, push the plug end of the connector cable fully into the output jack on the external neurostimulator.

Using the Therapy-stop key

The red Therapy-stop key can be used when you need to immediately turn off the external neurostimulator. You must use either the patient or clinician programmer to turn the external neurostimulator on again. The Therapy-stop key is not an on/off control.

- Press the red Therapy-stop key.

Disconnecting the connector cable from the external neurostimulator

1. Using the clinician programmer or patient programmer, turn the external neurostimulator off.
2. Disconnect the connector cable from the external neurostimulator.

Device care and storage

- Keep new batteries available.
- Use the patient programmer daily to check the external neurostimulator batteries.
- Replace low or depleted batteries.
- Handle the device and system components with care. Do not drop, strike or step on the device or system components.
- Do not dismantle or tamper with the device.
- Clean the outside of the device with a damp cloth when necessary. Mild household cleaners will not damage the device or labels.
- The battery contacts may be cleaned periodically with a cotton swab dampened with alcohol. Do not use a pencil eraser or sandpaper.
- Store the external neurostimulator at room temperature. Avoid extreme hot or cold temperatures and direct sunlight.
- The device and system components are not waterproof. Do not allow moisture to get inside the device or system components.
- Dispose of depleted batteries and worn out devices according to local requirements.

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The external neurostimulator contains no serviceable components. If the external neurostimulator requires repair or is nonfunctional, send it to the appropriate address.

USA

Medtronic, Inc.
Neuromodulation Division
Returned Product Analysis
MS N800
PO Box 1250
Minneapolis, MN 55440-9087

Filename Date Time
UC200xxxxx EN
4.625 x 6 inches (117 mm x 152 mm)

Medtronic Confidential
ImplantManual.xsl - IPGTemplate.fm
Template version 6.1: 01-09-2014

Europe, Africa, Middle East, and Asia-Pacific countries
Medtronic EOC
Medical Equipment Service Europe
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Tel. 31-455664880
Fax 31-455668028

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Contacts for specific countries are listed inside this cover.

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


M929110A074 Rev A

Printing instructions: Refer to doc#
A88898 for Printing Instructions

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










M929110A074 Rev A


Table of contents

Description	9
Package contents	9
Patient identification card	9
Device specifications	9
Declaration of conformity	12

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The Indications sheet provides information about indications and related information. For customers in Japan, the appropriate package insert provides information about indications.

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[USA] The clinical summary provides information about the clinical study results for the neurostimulation system.

Filename Date Time
UC20000000x EN
4.6 x 6 inches (117 mm x 152 mm)

Medtronic Confidential
ImplantManual.xsl - IPGTemplate.fm
Template version 6.3: 10-10-2019

Instructions for use 13

- Verifying neurostimulator operation 13
- Creating a pocket for the neurostimulator 13
- Connecting the extension to the neurostimulator 13
- Implanting the neurostimulator 15
- Checking system integrity 16
- Completing the implant procedure 16

8 English 37601 2020-01-01

2020-01-01

M929110A074 Rev A

Description

The Medtronic Activa™ PC Model 37601 Neurostimulator is part of a neurostimulation system for deep brain stimulation.

Package contents

- Neurostimulator
- Torque wrench
- Product literature
- **[USA]** Warranty card
- Registration form
- Patient identification card

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[USA] The patient identification card packaged with the device is temporary; a permanent card will be mailed to the patient when Medtronic receives the registration form.

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Device specifications

The neurostimulator is a multi-programmable device that delivers stimulation through 1 or 2 leads. The stimulation settings are stored in programs. A program is a specific combination of pulse width, rate, and amplitude settings acting on a specific electrode combination. Up to four programs can be combined into a group. When using more than one program, the pulses are delivered sequentially—first a pulse from one program, then a pulse from the next program.

Pulse width, amplitude, and electrode polarity for each program within the group can have different values. Rate, rate limits, SoftStart/Stop™ and cycling for each program within the group have the same values.

Table 1. Operating values for the Active PC Model 37601 neurostimulator

Programmable parameter	Operating range and resolution ^a
Number of defined groups	1 to 4
Number of programs per group	1 to 4
Amplitude (voltage mode)	0 to 10.5 V with 0.05-V or 0.1-V resolution
Amplitude (current mode)	0 to 25.5 mA with 0.1-mA resolution
Amplitude – upper patient limit	Tracking limit (by hemisphere): +0 to +2 (0.2 resolution); +2 to +4.5 (0.5 resolution)
Amplitude – lower patient limit	Tracking limit (by hemisphere): -0 to -2 (0.2 resolution); -2 to -4.0 (0.5 resolution); full range ^b
Pulse width	80 to 450 μ s (10- μ s resolution)
Pulse width – upper patient limit	Tracking limit: +0 to +100 μ s (10- μ s resolution)
Pulse width – lower patient limit	Tracking limit: -0 to -100 μ s (10- μ s resolution)
Rate (voltage mode)	2 to 250 Hz (resolution: 1 Hz from 2 Hz to 10 Hz, 5 Hz from 10 Hz to 250 Hz) ^c
Rate (current mode)	30 to 250 Hz (5-Hz resolution) ^c
Rate – upper patient limit	Tracking limit: +0 to +50 Hz (10-Hz resolution)
Rate – lower patient limit	Tracking limit: -0 to -50 Hz (10-Hz resolution)
SoftStart/Stop	Off, On: 1-, 2-, 4-, or 8-second ramp duration
Cycling	Off, On: 0.1 s to 24 hr (resolution: 0.1 s from 0.1 s to 1 s, 1 s from 1 s to 59 s, 1 min from 1 min to 59 min, 1 hr from 1 hr to 24 hr)

^a Interlocks will prevent the use of some parameter combinations.

^b Full range = -10.5 V (voltage mode); -25.5 mA (current mode).

^c Rate limited to 125 Hz when two programs are active on a single lead.

Table 2. Physical characteristics of the Active PC Model 37601 neurostimulator^a

Description	Value
Connector type	Octapolar, in-line 2.8-mm (0.110-in) spacing
Height	65.0 mm (2.6 in)
Length	49.0 mm (1.9 in)
Thickness	
case	15.0 mm (0.6 in)
connector	15.0 mm (0.6 in)
Weight	67.0 g (2.4 oz)
Volume	39.0 cm ³
Power source	6.3 Amp hours, 3.2 V HCSVO ^b cell
Temperature limitation	-18 to +52 °C (0 to +126 °F)
Serial number model designator ^c	NKM
Radiopaque Identification (ID) code	NKD
Transmitter	
Carrier frequency	175 kHz
Output level (at 300 m)	-45.8 dBuV/m

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^b Hybrid combined silver vanadium oxide.

^c The serial number is the model designator followed by a number. The clinician programmer displays the entire serial number beginning with the model designator.

Filename Date Time
UC200XXXXX EN
4.6 x 6 inches (117 mm x 152 mm)

Medtronic Confidential
ImplantManual.xsl - IPGTemplate.fm
Template version 8.3: 10-10-2019

Table 3. Material of components in the Activa PC Model 37601 package

Components	Material	Material contacts human tissue
Neurostimulator		
Case	Titanium	Yes
Connector block	Polyurethane, silicone rubber, silicone medical adhesive	Yes
Grommets, seals	Silicone rubber	Yes
Setscrews	Titanium alloy	Yes
Adhesive	Silicone medical adhesive	Yes

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M929110A074 Rev A

Instructions for use

Implanting physicians should be experienced in stereotactic and functional neurosurgery and have expertise with functional stereotactic neurosurgical treatment of movement disorders and deep brain stimulation procedures, as well as thoroughly familiar with all product labeling.

△ Cautions:

- When using sharp instruments near the neurostimulator, be extremely careful to avoid nicking or damaging the case or the connector block. Damaging the neurostimulator may require surgical replacement.
- Do not use saline or other ionic fluids at connections, which could result in a short circuit.

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Before opening the sterile neurostimulator package, verify that the neurostimulator is operable by using the clinician programmer to interrogate the neurostimulator and read the neurostimulator battery service life. (Refer to the software manual for instructions on how to read the battery service life.)

△ **Caution:** Do not implant a neurostimulator if it was dropped onto a hard surface from a height of 30 cm (12 in) or more, because the neurostimulator may be damaged and fail to operate properly.

Note: The neurostimulator pocket may be flushed with a nonionic antibiotic solution; do not submerge the neurostimulator in fluid.

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△ **Caution:** Before connecting components, wipe off any body fluids and dry all connections. Fluids in the connection may result in stimulation at the connection site, intermittent stimulation, or loss of stimulation.

1. Wipe the extension connectors with sterile gauze. If necessary, use sterile (United States Pharmacopeia [USP]) water or a nonionic antibiotic solution.
2. Make sure the connector block receptacles are dry and clean.
3. If using two leads, identify the left and right lead-extensions.

Note: It may be necessary to tug gently on either lead-extension and view its movement at the lead site to identify left and right lead-extensions.

4. Insert each extension connector into the appropriate neurostimulator socket until it is fully seated within the connector block (Figure 1).

Notes:

- Insert the identified left extension into the front (0-7) neurostimulator socket (socket I) and the right extension into the back (8-15) neurostimulator socket (socket II) to simplify electrode assignment during programming.
- During insertion, some resistance is typical.
- To retract the setscrews, insert the torque wrench into the self-sealing grommet and rotate the setscrews counterclockwise; however, do not remove the setscrews from the connector block.

⚠ **Caution:** Do not insert the extension connector into the connector block if the setscrews are not sufficiently retracted. If the setscrews are not retracted, the setscrews may damage the extension and the extension will not be seated fully into the connector block.


Figure 1. Insert the extension connectors fully into the neurostimulator socket.

Note: Insert a neurostimulator plug (from the extension kit) into an unused neurostimulator socket.

5. For each extension or plug, fully insert the torque wrench (packaged with the neurostimulator) into each self-sealing grommet of the connector block and tighten each setscrew (Figure 2).

△ Cautions:

- To prevent undertightening the neurostimulator setscrews, do not use the torque wrench from the extension kit. Undertightening may result in insufficient electrical contact within the connector block, which may cause intermittent stimulation.
- Be sure the torque wrench is fully inserted into the self-sealing grommet. If the torque wrench is not fully inserted, the grommet may be damaged, resulting in inadequate stimulation.
- Before tightening setscrews, ensure that the extension connector is fully inserted into the connector block to prevent damaging the extension.
- Verify that each leaf of the self-sealing grommet is closed after the torque wrench is withdrawn. If fluid leaks through a grommet seal that is not fully closed, the patient may experience shocking, burning, or irritation at the neurostimulator implant location, or intermittent stimulation, or loss of stimulation.


Figure 2. Tightening the setscrews in the self-sealing grommet.

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△ Cautions:

- Ensure that the neurostimulator is placed no deeper than 4 cm (1.5 in) below the skin and is parallel to the skin. If the neurostimulator is too deep or is not parallel to the skin, telemetry may be unsuccessful.
- Do not coil excess extension in front of the neurostimulator. Wrap excess extension around the perimeter (Figure 3) of the neurostimulator to minimize subcutaneous pocket depth, help minimize potential damage during

neurostimulator replacement surgery, help minimize potential kinking of the extension, and minimize interference with telemetry. Excess extension should not exceed two wraps around the perimeter of the neurostimulator. Extension lengths requiring more than two wraps can interfere with telemetry.


Figure 3. Wrap excess extension around the perimeter of the neurostimulator.

2. Use the suture holes in the connector block to secure the neurostimulator to the muscle fascia with nonabsorbable silk.

Checking system integrity

Warning: To use the programming components (external and nonsterile) in a sterile field, place a sterile barrier between the patient and the programming components to prevent infection. Do not sterilize any of the programming components. Sterilization may damage the programming components.

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1. To ensure proper connection of each extension to the neurostimulator, use the clinician programmer to program left and right lead assignments and the basic stimulation parameters, check the battery status, and check the electrode impedances to rule out a short or open circuit.
2. If the system integrity test results are not acceptable, refer to "Connecting the extension to the neurostimulator" on page 13.

Completing the implant procedure

1. Close and dress all incisions.
2. Ensure that a patient control device is given to the patient.

Filename Date Time
UG200xxxxx EN
4.6 x 6 inches (117 mm x 152 mm)

Medtronic - Confidential
ImplantManual.xsl - IPGTemplate.fm
Template version 6.3: 10-10-2019

3. Complete the device tracking and patient registration paperwork and return the documents to Medtronic.

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Authorized Representative 
In the European Community

Medtronic B.V.
Earl Bakkenstraat 10,
6422 PJ Heerlen,
The Netherlands
Tel. +31-45-586-8000

Europe/Africa/Middle East Headquarters

Medtronic International Trading Sàrl
Route du Molliu 31,
Case Postale 84
CH - 1131 Tolochenaz,
Switzerland
www.medtronic.eu
Tel. +41-21-802-7000

Asia-Pacific
Medtronic International Ltd.
50 Pasir Panjang Road,
#04-51 Mapletree Business City,
Singapore 117384,
Singapore
Tel. +65-6870-5510

© Medtronic 2019
All Rights Reserved

2020-01-01
M929110A074 Rev A

M929110A074 Rev A

2020-01-01

Filename Date Time
UC200cccc EN
4.6 x 6 inches (117 mm x 152 mm)

Medtronic Confidential
ImplantManual.usl - IPGTemplate.fm
Template version 6.3: 10-10-2019

Medtronic

Activa™ RC
Multi-program rechargeable
neurostimulator

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
CE 0123
2008

M927170A071 Rev A

Printing instructions: Refer to doc#
A88896 for Printing Instructions



2020-01-01

Filename Date Time
UC200xxxxx EN
4,6 x 6 inches (117 mm x 152 mm)


Medtronic Confidential
ImplantManual.xsl - IPGTemplate.fm
Template version 6.3: 10-10-2019

M827170A071 Rev A

Printing instructions: Refer to doc#
A88896 for Printing Instructions

2020-01-01

Explanation of symbols on product or package labeling

Refer to the appropriate product for symbols that apply.

Open here

Do not reuse

Sterilized using ethylene oxide

Date of manufacture

Manufacturer

Use by

Serial number

Consult instructions for use

Conformité Européenne (European Conformity). This symbol means that the device fully complies with AIMD Directive 90/385/EEC (NB 0123).

Authorized representative in the European Community

Do not resterilize

Do not use if package is damaged



For USA audiences only

PIN number

Magnetic Resonance (MR) Conditional

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

M927170A071 Rev A

Table of contents

Description	7
Package contents	7
Patient identification card	7
Device specifications	7
Declaration of conformity	10

Information available for the Deep Brain Stimulation Therapies:

The information for prescribers manual provides information about contraindications, warnings, precautions, adverse events, sterilization, and component disposal. The information for prescribers addendums contain indication-specific information. For customers in Japan, the appropriate package insert provides information about safety, contraindications, warnings, precautions, and adverse events.

The indications sheet provides information about indications and related information. For customers in Japan, the appropriate package insert provides information about indications.

The system eligibility and battery longevity manual describes programming considerations and provides battery longevity information to aid in the appropriate neurostimulator selection.

MRI guidelines provide information about any MRI conditions and MRI-specific contraindications, warnings, and precautions for MRI scans with the neurostimulation system.

Product manuals, such as programming guides, recharging guides, and implant manuals provide device descriptions, package contents, device specifications, product-specific warnings and precautions, and instructions for use.

USA The clinical summary provides information about the clinical study results for the neurostimulation system.

Filename Date Time
UC200xxxxx EN
4.6 x 6 inches (117 mm x 152 mm)

Medtronic Confidential
ImplantManual.xsl - IPGTemplate.fm
Template version 6.3: 10-10-2019

Instructions for use 11

- Charging the neurostimulator battery 11
- Verifying neurostimulator operation 11
- Creating a pocket for the neurostimulator 12
- Connecting the extension to the neurostimulator 12
- Implanting the neurostimulator 14
- Checking system integrity 15
- Completing the implant procedure 15

6 English 37612 2020-01-01

2020-01-01

M927170A071 Rev A

Description

The Medtronic Activa™ RC Model 37612 Neurostimulator is part of a neurostimulation system for deep brain stimulation.

Package contents

- Neurostimulator
- Torque wrench
- Product literature
- Warranty card
- Registration form
- Patient identification card

Patient identification card

A patient identification card is packaged with this device. Advise the patient to carry the identification card at all times.

USA: The patient identification card packaged with the device is temporary; a permanent card will be mailed to the patient when Medtronic receives the registration form.

The implant registration form registers the device warranties and creates a record of the device in Medtronic's implant data system.

Device specifications

The neurostimulator is a multi-programmable, rechargeable device that delivers stimulation through 1 or 2 leads. The stimulation settings are stored in programs. A program is a specific combination of pulse width, rate, and amplitude settings acting on a specific electrode combination. Up to four programs can be combined into a group, with a maximum of 2 programs per lead. When using more than one program, the pulses are delivered sequentially—first a pulse from one program, then a pulse from the next program.

Pulse width, amplitude, and electrode polarity for each program within the group can have different values. Rate, rate limits, ramping, and cycling for each program within the group have the same values.

Table 1. Operating values for the Activa RC Model 37612 neurostimulator

Programmable parameter	Operating range and resolution ^a
Number of defined groups	1 to 4 ^b
Number of programs per group	1 to 4 ^b
Amplitude (voltage mode)	0 to 10.5 V with 0.05 V or 0.1-V resolution
Amplitude (current mode)	0 to 25.5 mA with 0.1-mA resolution
Amplitude – upper patient limit	By hemisphere: Tracking limit: +0 to +2 (0.2 resolution); +2 to +4.5 (0.5 resolution)
Amplitude – lower patient limit	By hemisphere: Tracking limit: -0 to -2 (0.2 resolution); -2 to -4.0 (0.5 resolution); full range ^c
Pulse width	60 to 450 μ s (10- μ s resolution)
Pulse width – upper patient limit	Tracking limit: +0 to +100 μ s (10- μ s resolution)
Pulse width – lower patient limit	Tracking limit: -0 to -100 μ s (10- μ s resolution)
Rate (voltage mode)	2 to 250 Hz (resolution: 1 Hz from 2 Hz to 10 Hz, 5 Hz from 10 Hz to 250 Hz) ^d
Rate (current mode)	30 to 250 Hz (5-Hz resolution) ^d
Rate – upper patient limit	Tracking limit: +0 to +50 Hz (10-Hz resolution)
Rate – lower patient limit	Tracking limit: -0 to -50 Hz (10-Hz resolution)
SoftStart/Stop™	Off, On: 1-, 2-, 4-, or 8-second ramp duration
Cycling	Off, On: 0.1 s to 24 hr (resolution: 0.1 s from 0.1 s to 1 s, 1 s from 1 s to 59 s, 1 min from 1 min to 59 min, 1 hr from 1 hr to 24 hr)

^a Interlocks will prevent the use of some parameter combinations.

^b No more than 16 programs may be defined within the 4 groups.

^c Full range = -10.5 V (voltage mode); -25.5 mA (current mode)

^d Rate limited to 125 Hz when two programs are active on a single lead.

Table 2. Physical characteristics of the Activa RC Model 37612 neurostimulator^a

Description	Value
Connector type	Octapolar, in-line 2.8-mm (0.110-in) spacing
Height	54 mm (2.1 in)
Length	54 mm (2.1 in)
Thickness	
case	9 mm (0.4 in)
connector	11 mm (0.4 in)
Weight	40 g (1.4 oz)
Volume	22 cm ³
Battery life	15 years (9 years as shipped) ^b
Power source ^c	Lithium ion rechargeable battery
Serial number model designator ^d	NKG
Radiopaque Identification (ID) code	NKG
Transmitter	
Carrier frequency	175 kHz
Output level (at 300 m)	-48 dBuV/m

^a All measurements are approximate.

^b Service life extension is required before 9-year EOS.

^c The neurostimulator is not shipped with a full battery charge. Refer to the charging system user manual for neurostimulator charging instructions.

^d The serial number is the model designator followed by a number. The clinician programmer displays the entire serial number beginning with the model designator.

Filename Date Time
UC200xxxxx EN
4.6 x 6 inches (117 mm x 152 mm)

Medtronic Confidential
ImplantManual.xsl - JPGTemplate.fm
Template version 6.3: 10-16-2019

Table 3. Material of components in the Active RC Model 37612 package

Components	Material	Material contacts human tissue
Neurostimulator		
Case	Titanium	Yes
Connector block	Polysulfone, silicone rubber, silicone medical adhesive	Yes
Grommets, seals	Silicone rubber	Yes
Screws	Titanium alloy	Yes
Adhesive	Silicone medical adhesive	Yes

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For additional information, contact the appropriate Medtronic representative listed on the back cover of this manual.

10 English 37612 2020-01-01

2020-01-01

M927170A071 Rev A

Instructions for use

Implanting physicians should be experienced in stereotactic and functional neurosurgery and have expertise with functional stereotactic neurosurgical treatment of movement disorders and should be thoroughly familiar with all product labeling.

⚠ **Warning:** DO NOT use the recharger on an unhealed wound. The recharging system is not sterile, and contact with the wound can cause an infection.

⚠ **Warning:** Do not use saline or other ionic fluids at connections, which could result in a short circuit.

⚠ **Caution:** Advise patients to charge the neurostimulator on a regular basis, taking into consideration the therapeutic parameters and battery usage of the patient, to prevent the battery from overdischarging. If the neurostimulator battery is allowed to overdischarge, the patient cannot charge the neurostimulator; however, the clinician may be able to restore the battery function using the Physician Recharge Mode on the recharger (refer to the troubleshooting section of the software manual).

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- Battery function is restored; however, charging sessions may be more frequent because battery capacity has been reduced.
- Battery function is not restored and the neurostimulator must be surgically replaced. Battery function is not restored when:
 - the neurostimulator battery is permanently damaged.
 - the neurostimulator battery has been overdischarged and restored twice before. The third time the battery is overdischarged, the neurostimulator will reach end of service. Surgery is required to replace the neurostimulator.

⚠ **Caution:** When using sharp instruments near the neurostimulator, be extremely careful to avoid nicking or damaging the case or the connector block. Damaging the neurostimulator may require surgical replacement.

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Charge the neurostimulator battery before opening the package. For charging instructions, refer to the charging system user manual.

Verifying neurostimulator operation

Before opening the sterile neurostimulator package, verify that the neurostimulator is operable by using the clinician programmer to interrogate the neurostimulator and read the neurostimulator battery charge level. (Refer to the software manual for instructions on how to read the battery charge level.)

⚠ **Caution:** Do not implant a neurostimulator if it was dropped onto a hard surface from a height of 30 cm (12 in) or more, because the neurostimulator may be damaged and fail to operate properly.

Note: The neurostimulator pocket may be flushed with a nonionic antibiotic solution.

⚠ Caution: Do not submerge the neurostimulator in fluid. Fluids in the connection may result in stimulation at the connection site, intermittent stimulation, or loss of stimulation.

Creating a pocket for the neurostimulator

1. Make an incision slightly larger than the width of the neurostimulator.
2. Create a subcutaneous pocket for the neurostimulator.
3. Verify the pocket size.

Connecting the extension to the neurostimulator

⚠ Caution: Before connecting components, wipe off any body fluids and dry all connections. Fluids in the connection may result in stimulation at the connection site, intermittent stimulation, or loss of stimulation.

1. Wipe the extension connectors with sterile gauze. If necessary, use sterile (United States Pharmacopeia [USP]) water or a nonionic antibiotic solution.
2. Make sure the connector block receptacles are dry and clean.
3. Insert each extension connector into the appropriate neurostimulator socket until it is seated fully within the connector block (Figure 1).

Notes:

- During insertion, some resistance is typical.
- To retract the setscrews, insert the torque wrench into the self-sealing grommet and rotate the setscrews counterclockwise; however, do not remove the setscrews from the connector block.

⚠ Warning: Do not insert the extension connector into the connector block if the setscrews are not sufficiently retracted. If the setscrews are not retracted, the setscrews may damage the extension and the extension will not be seated fully into the connector block.


Figure 1. Insert the extension connectors into the neurostimulator.

Note: Insert a connector plug (from an accessory kit) into unused neurostimulator socket.

4. For each extension or plug, fully insert the torque wrench (packaged with the neurostimulator) into each self-sealing grommet of the connector block and tighten each setscrew (Figure 2).

Warnings:

- Be sure the torque wrench is fully inserted into the self-sealing grommet. If the torque wrench is not fully inserted, the grommet may be damaged, resulting in inadequate stimulation.
- Before tightening setscrews, ensure that the extension connector is fully inserted into the connector block to prevent damaging the extension.

Cautions:

- To prevent undertightening the neurostimulator setscrews, do not use the torque wrench from the extension kit. Undertightening may result in insufficient electrical contact within the connector block, which may cause intermittent stimulation.
- Verify that each leaf of the self-sealing grommet is closed after the torque wrench is withdrawn. If fluid leaks through a grommet seal that is not fully closed, the patient may experience shocking, burning, or irritation at the neurostimulator implant location, or intermittent stimulation, or loss of stimulation.


Figure 2. Tightening the setscrews in the self-sealing grammet.

Implanting the neurostimulator

1. Place the neurostimulator into the subcutaneous pocket with the Medtronic logo facing outward, away from muscle tissue, and ensure that the extension is not bent sharply.

Warnings:

- Ensure that the neurostimulator is placed no deeper than 1 cm (0.4 in) below the skin and is parallel to the skin. If the neurostimulator is too deep or is not parallel to the skin, recharge may be inefficient or unsuccessful.
- Do not coil excess extension in front of the neurostimulator. Wrap excess extension around the perimeter (Figure 3) of the neurostimulator to minimize subcutaneous pocket depth, help minimize potential damage during neurostimulator replacement surgery, help minimize potential kinking of the extension, and minimize interference with telemetry and recharge operation. Excess extension should not exceed two wraps around the perimeter of the neurostimulator. Extension lengths requiring more than two wraps can interfere with telemetry.


Figure 3. Wrap excess extensions around the perimeter of the neurostimulator.

⚠ **Caution:** Position the neurostimulator with the Medtronic logo facing outward. If implanted with the Medtronic logo facing inward, the neurostimulator will be difficult to change.

2. Secure the neurostimulator to the muscle fascia with nonabsorbable silk.

⚠ **Warning:** Use the suture holes in the connector block to secure the neurostimulator to the muscle fascia with nonabsorbable silk. Device inversion may damage the extension, which may cause inadequate stimulation and require a surgical revision or replacement.

Checking system integrity

⚠ **Warning:** To use the programming components (external and nonsterile) in a sterile field, place a sterile barrier between the patient and the programming components to prevent infection. Do not sterilize any of the programming components. Sterilization may damage the programming components.

Note: The neurostimulator should be in the pocket during system interrogation for integrity to ensure proper readings.

1. To ensure proper connection of each extension to the neurostimulator, use the clinician programmer to program the basic stimulation parameters, check the battery status, and check the electrode impedances to rule out a short or open circuit.
2. If the system integrity test results are not acceptable, refer to "Connecting the extension to the neurostimulator" on page 12.

Completing the implant procedure

1. Close and dress all incisions.

Filename Date Time
UC200xxxxx EN
4.6 x 6 inches (117 mm x 152 mm)

Medtronic Confidential
ImplantManual.xsl - IPGTemplate.fm
Template version 6.3: 10-10-2019

2. Ensure that a patient control device is given to the patient.
3. Complete the device tracking and patient registration paperwork and return the documents to Medtronic.

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
M992991A008 Rev A 2019-04-01



Printing instructions. Refer to doc# A88898 for Printing
Instructions

Filename Date Time
UC200xxxxxx EN
4.6 x 6 inches (117 mm x 152 mm)


Medtronic Confidential
EMCDeclarations.xsl - EMCDeclarationsTemplate.fm
Template version 5.2: 02-07-2019

M692991A008 Rev A 2019-04-01

Printing instructions: Refer to doc# A88896 for Printing
Instructions

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



4 English 37642 2019-04-01

M992991A006 Rev A 2019-04-01

Wireless quality of service (QoS) and operating characteristics

Near-field magnetic induction communication

The Medtronic Model 37642 patient programmer uses near-field magnetic induction communication to interact with compatible medical devices.

Hold the Model 37642 patient programmer directly over the neurostimulator with the screen facing outward while adjusting therapy.

Alternatively, the Model 37092 detachable antenna can be used for viewing the patient programmer screen while adjusting therapy. Hold the antenna or attach it to clothing directly over the neurostimulator.

If the wireless connection is lost during a programming session, the user should reposition the programmer or the detachable antenna over the neurostimulator to reestablish the connection.

See Table 1 for the wireless operating characteristics.

Table 1. Wireless operating characteristics: Near-Field Magnetic Induction Communication for Model 37642

Radio type	Near-field magnetic induction communication
Operating frequency band	175 kHz
Effective radiated power output	H-field strength of <30 dBµA/m at 10 m
Receiving bandwidth	Up to 50 kHz (occupied bandwidth of the transmitter)
Modulation (modulation is model-specific)	On-off keying (OOK) burst
Effective range	With external neurostimulator (ENS): 2.5 cm (0.98 in) – approximate maximum distance between the Model 37642 patient programmer and the ENS. With implantable neurostimulator (INS): 1.5 cm (0.59 in) – approximate maximum distance of the Model 37642 patient programmer from the surface of the skin directly over the INS; distance may vary depending on the implant depth of the INS.

Wireless security

Near-field magnetic induction communication

The near-field magnetic induction communication uses a model-specific device access code to initiate a programming session. Hold the Model 37642 patient programmer or the detachable antenna directly over the neurostimulator. This close contact will mitigate inadvertent programming, malicious programming, and eavesdropping.

Electromagnetic compatibility (EMC) declaration

The Medtronic Model 37642 patient programmer (also referred to herein as "external device") is portable medical electrical equipment and has been tested and conforms to the requirements of IEC 60601-1-2:2014.

Medical electrical equipment needs special precautions regarding electromagnetic compatibility (EMC) and needs to be put into service according to the EMC information provided in this EMC declaration.

- For the environments in which the Model 37642 patient programmer is intended to be used, refer to the information in Table 2 on page 7 of this EMC declaration.
- The Model 37642 patient programmer and all accessories have no essential performance.
- **Note:** The following information does not apply to equipment that Medtronic identifies as compatible and intended to be used with the external device:
 - Portable and mobile radio-frequency (RF) communications equipment (for example, mobile or cellular phones and amateur radio equipment) can affect medical electrical equipment.
 - Portable radio-frequency (RF) communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Model 37642 patient programmer, including cables specified by Medtronic. Otherwise, degradation of the performance of this external device could result.
 - Use of this external device adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this external device and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers, and cables other than those specified or provided by Medtronic could result in increased electromagnetic emissions or decreased electromagnetic immunity of this external device and result in improper operation.
- Refer to the manual for the Model 37642 patient programmer for the medical devices, equipment, and accessories that Medtronic identifies as compatible with this external device.
- Other equipment could interfere with the external device, even if the other equipment complies with CISPR emission requirements.

- This external device is MR Unsafe and should not be allowed in the MRI scanner (magnet) room. This external device contains ferromagnetic material, which can be affected by the MRI magnet.
- The use of this external device can be affected by interference when in the vicinity of equipment that uses RF or microwave energy (such as diathermy, RF ablation, or an RFID reader).

For more information on MRI, diathermy, and other sources of electromagnetic interference (EMI), refer to other manuals provided with the patient's Medtronic neuromodulation system.

Tables 2, 3, and 4 apply to the Model 37642 patient programmer, which is an external device internally powered by nonrechargeable batteries.

Table 2. (IEC 60601-1-2:2014) Examples of environments of intended use

Electromagnetic environments	Examples
Home healthcare environment	Restaurants, cafes, shops, stores, markets, schools, churches, libraries, outdoors (streets, sidewalks, parks), domiciles (residences, homes, nursing homes), vehicles (cars, buses, trains, boats, planes, helicopters), train stations, bus stations, airports, hotels, hostels, pensions, museums, theatres
Professional healthcare facility environment	Physician offices, dental offices, clinics, limited care facilities, freestanding surgical centers, freestanding birthing centers, multiple treatment facilities, hospitals – emergency rooms, patient rooms, intensive care, outside the RF shielded room of a medical electrical (ME) system for magnetic resonance imaging (MRI), surgery rooms except near high-frequency (HF) surgical equipment

Table 3. Electromagnetic emissions: Test information and results

Emissions test method	IEC 60601-1-2: 2014 Test level requirement and compliance test level	Test applicability, guidance, and comments
		Model 37642 patient programmer (battery-powered, nonrechargeable)
RF conducted emissions CISPR 11 (IEC 55011)	Group 1, Class B	Not applicable: external device is internally powered by batteries and has no patient port.
RF radiated emissions CISPR 11 (IEC 55011)	Group 1, Class B	Applicable and compliant.
Harmonic emissions IEC 61000-3-2	Test levels are according to classification of equipment (Class A, B, C, D) as listed in clause 7 of IEC 61000-3-2.	Not applicable: external device is internally powered by batteries.
Voltage fluctuations/flicker emissions IEC 61000-3-3	Test levels are according to clause 5 of IEC 61000-3-3.	Not applicable: external device has no AC port or DC port.

8 English 37642 2019-04-01

M992991A008 Rev A 2019-04-01

Table 4. Electromagnetic immunity: Testing information and results

Immunity test method	IEC 60601-1-2: 2014 Test level requirement and compliance test level	Test applicability, guidance, and comments
		Model 37642 patient programmer (battery-powered, nonrechargeable)
Radiated RF IEC 61000-4-3	10 V/m: 80 MHz to 2.7 GHz 80% AM at 1 kHz	Applicable and compliant.
Radiated RF – Proximity to portable and mobile radio communication equipment IEC 61000-4-3	Exposure to 385 MHz to 5785 MHz 385 MHz: 27 V/m 450 MHz: 28 V/m 710, 745, 780 MHz: 9 V/m 810, 870, 930 MHz: 28 V/m 1720, 1845, 1970 MHz: 28 V/m 2450 MHz: 28 V/m 5240, 5500, 5800 MHz: 9 V/m	Applicable and compliant. Recommended minimum separation distance is 30 cm (12 inches).
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	Applicable and compliant.
Conducted disturbances induced by RF IEC 61000-4-6	3 V rms: 150 kHz to 80 MHz, 80% AM at 1 kHz 6 V rms in ISM and amateur radio bands from 150 kHz to 80 MHz	Not applicable: external device has no patient port, no AC port, no SIP/SOP port.

Table 4. Electromagnetic immunity: Testing information and results (continued)

Immunity test method	IEC 60601-1-2: 2014 Test level requirement and compliance test level	Test applicability, guidance, and comments
		Model 37642 patient programmer (battery-powered, nonrechargeable)
Electrostatic discharge (ESD) IEC 61000-4-2	• Contact: ± 8 kV • Air: ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV The following are all the test levels that were tested: • Contact: ± 2 kV, ± 4 kV, ± 8 kV • Air: ± 2 kV, ± 4 kV, ± 8 kV, ± 10 kV, ± 12 kV, ± 15 kV	Applicable and compliant. Enclosure and all ports were evaluated.
Electrical fast transient/burst IEC 61000-4-4	• ± 2 kV for power supply lines • ± 1 kV for input/output lines	Not applicable: external device has no patient port, no AC port, no SIP/SOP port.
Surge IEC 61000-4-5	• ± 0.5 kV, ± 1 kV differential mode • ± 0.5 kV, ± 1 kV, ± 2 kV common mode	Not applicable: external device has no patient port, no AC port, no SIP/SOP port.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	• 0% U_T for 0.5 cycle • 0% U_T for 1 cycle • 70% U_T for 25/30 cycles • 0% U_T for 250/300 cycles NOTE: U_T is AC mains voltage prior to application of test level.	Not applicable: external device has no AC port.

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



Authorized Representative 
In the European Community
Medtronic B.V.
Earl Bakkenstraat 10,
6422 PJ Heerlen,
The Netherlands
Tel. +31-45-566-8000

Europe/Africa/Middle East Headquarters
Medtronic International Trading Sàrl
Route du Molliat 31,
Case Postale 84
CH - 1131 Tolochenaz,
Switzerland
www.medtronic.eu
Tel. +41-21-802-7000

Asia-Pacific
Medtronic International Ltd.
Suits 1106-11, 11/F, Tower 1, The Gateway,
25 Canton Road, Tsimshatsui,
Kowloon,
Hong Kong
Tel. +852-2919-1300

© Medtronic, Inc. 2019
All Rights Reserved
2019-04-01
M992991A008 Rev A

M992991A008 Rev A 2019-04-01

Filename Date Time
UC200xxxxx EN
4.6 x 6 inches (118 mm x 152 mm)

Medtronic Confidential
PPManual.xls! - PatientProgrammerTemplate.fm
Template version 5.3: 09-13-2019

Medtronic

PATIENT PROGRAMMER

37642

Medtronic DBS Therapy user manual
Activa™ PC Model 37601
Activa™ RC Model 37612
Activa™ SC Model 37602
Activa™ SC Model 37603

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CE 0123
2008

M929532A059 Rev A 2020-01-01

Printing instructions: Refer to doc# A88898 for
Printing Instructions

Filename Date Time
UC200xxxxx EN
4.6 x 6 inches (116 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.3: 09-13-2019

Medtronic and Medtronic logo are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.

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M029532A050 Rev A 2020-01-01

Label symbols

Explanation of symbols on products and packaging. Refer to the appropriate product to see symbols that apply.

CE 0123 Conformité Européenne (European Conformity). This symbol means that the device fully complies with AIMD Directive 90/385/EEC (NB 0123).

Consult instructions for use

Manufacturer

Temperature limitation

Keep dry

Serial number

PIN number

Authorized representative in the European community

For USA audiences only

IEC 60601-1/EN60601-1, Type BF Equipment

Label symbols

37642 2020-01-01 English

3

M929532A059 Rev A 2020-01-01

Antenna jack

Magnetic Resonance (MR)
Conditional

Magnetic Resonance (MR) Unsafe

Do not dispose of this product in the unsorted municipal waste stream. Dispose of this product according to local regulations. See <http://recycling.medtronic.com> for instructions on proper disposal of this product.

Chinese Standard (SJ/T11364-2006)
Logo: Electronic Information Products Pollution Control Symbol. (The "e" in this logo signifies none of the restricted substances are present above permitted levels.)

The RCM Certification mark signifies that this device complies with Australia/New Zealand EMC requirements for electrical and electronic equipment.

Label symbols

4 English 37642 2020-01-01

M929532A059 Rev A 2020-01-01

Table of contents

Label Symbols 3

Glossary 10

1 Introduction 18

How to use this manual 18

For important safety information 20

Purpose of your patient programmer 20

Purpose of the neurostimulation system
(indications) 21

Precautions 21

2 Using your patient programmer 26

How the patient programmer works 26

Multiple neurostimulators 29

Simple and advanced therapy
modes 29

Summary of keys 32

Using the Navigator key 36

The Therapy screen 38

Table of contents

37642 2020-01-01 English 5

Synchronizing the patient programmer
and your neurostimulator 41

Turning your therapy on or off 44

The neurostimulator battery 46

About rechargeable batteries 46

Checking the nonrechargeable battery
status 47

Possible nonrechargeable battery status
conditions 48

Customizing your patient
programmer 50

Changing preferences 51

Patient programmer alert 57

Patient programmer accessories 58

Carrying case and patient programmer ID
label 58

Detachable antenna 60

Using the detachable antenna 61

3 Adjusting therapy settings with Advanced mode 66

Adjusting your therapy settings 66

Therapy settings 68

Changing therapy settings in Advanced
mode 69

Selecting a new group 73
Returning your therapy to original
settings 77
Understanding therapy parameter
limits 81

4 MRI examinations 84

If you have an MRI appointment 84
Responsibilities of the patient in preparing
for the MRI appointment 84
Preparing your neurostimulation system
for the MRI scan 85
Turning off therapy before the MRI
scan 86
 Using the patient programmer to turn off
 therapy 87
 Using the clinician programmer to turn off
 therapy 88
Turning therapy back on after the MRI
scan 89
Activating settings to allow therapy on
during the MRI scan 90
 Using the patient programmer to activate
 a new group 91

Table of contents

37642 2020-01-01 English 7

Using the clinician programmer to activate
a new group 92

Returning therapy to your original group
setting after the MRI scan 94

5 Important information about your rechargeable neurostimulator 98

The rechargeable neurostimulator
battery 98

Check and charge: make it a habit 100

Checking the rechargeable
neurostimulator battery 102

Low battery status 103

Rechargeable battery overdischarge 105

Charging the neurostimulator battery 106

Patient programmer alert 108

6 Troubleshooting 110

Patient Programmer screens 110

Warning screens 110

Communication screen 113

Information screens 114

Table of contents

8 English 37642 2020-01-01

M929532A059 Rev A 2020-01-01

Possible problems and solutions 121

7 Maintenance 128

Patient programmer batteries 128

Checking the patient programmer
batteries 129

Replacing the patient programmer
batteries 132

Cleaning and care 134

Safety and technical checks 135

Battery and patient programmer
disposal 136

Declaration of conformity 136

Specifications 136

8 Assistance for the patient programmer 140

Index 142

Table of contents

37642 2020-01-01 English 9

Glossary

Active Group - Available in Advanced mode only, the active group is the current group selected for your neurostimulator. See "Group" for more information.

Amplitude - The strength or intensity of an electrical pulse.

Caution - A statement describing actions that could result in damage to or improper functioning of a device.

Charging system - Equipment used to charge the battery inside a rechargeable implanted neurostimulator.

Clinician - A healthcare professional such as a doctor or nurse.

Clinician programmer - A device used by a clinician to send instructions to a neurostimulator.

Contraindication - A condition or circumstance when a person should not have a neurostimulation system.

Glossary

10 English 37642 2020-01-01

M929532A059 Rev A 2020-01-01

Cycling - See Stimulation cycling.

Deep Brain Stimulation (DBS) - The delivery of electrical pulses to a targeted area or areas of the brain.

Discharged battery status (rechargeable battery only) - The rechargeable battery is depleted and should be charged as soon as possible. When the battery is in a discharged state, therapy is not available.

Electrode - A metal piece near the tip of the lead. Electrodes deliver electrical pulses to a targeted area or areas of the brain.

Electromagnetic interference (EMI) - A strong field of energy near electrical or magnetic devices that could prevent the neurostimulator from functioning properly.

EOS (End of Service) - A notification that the neurostimulator has reached its end of service. At EOS, the neurostimulator no longer delivers the electrical pulses that provide therapy.

ERI (Elective replacement indicator) - A notification that the neurostimulator is nearing its end of service.

Group - A group is a collection of therapy settings on your patient programmer that your doctor creates for you. Each group may be defined for a different activity, symptom, or time of day. Available in Advanced mode only, this feature allows you to adjust your therapy for different situations or for personal preference.

Group status row - The Parameter row on the Therapy screen displays the active group.

Implantable neurostimulator (INS) - See Neurostimulator.

Indication - The purpose of the neurostimulation system and the medical condition for which it may be implanted.

Lead - The lead is a set of thin wires covered with a protective coating. It carries the therapy signal to the electrodes that deliver stimulation to the brain tissue.

Magnetic Resonance Imaging (MRI) - A type of medical procedure that scans your body using magnetic fields to provide detailed pictures of your anatomy.

Neurostimulation system - Components that deliver, control, and maintain electrical pulses to provide therapy to the brain.

Neurostimulator - The neurostimulator is the implanted device that generates and controls your DBS therapy.

Overdischarge (rechargeable battery only)
- The rechargeable battery cannot be recharged when it is overdischarged. Once the battery is in a discharged state, if it is not recharged it continues to lose charge. Eventually, the battery loses enough charge to permanently affect the neurostimulator. If this occurs, the battery is overdischarged.

Parameter - One of three therapy settings that adjust the electrical pulse: amplitude, pulse width, and rate of the electrical stimulation.

Parameter/Group row - The bottom row on the Therapy screen showing the parameter and active group settings (only available in Advanced mode).

Parameter settings - See Therapy settings.

Patient programmer - A hand-held device that allows you to turn your therapy on and off. It may also be used to adjust some therapy settings.

Precaution - See Caution.

Preferences - Adjustable settings on your patient programmer including audio, contrast, and text or icon display format. Preferences do not affect your therapy.

Pulse width - The length or duration of an electrical pulse.

Rate - The number of electrical pulses delivered each second.

Glossary

14 English 37642 2020-01-01

M929532A059 Rev A 2020-01-01

SoftStart/Stop - This feature, programmed by your clinician, starts and stops your therapy gradually by slowly increasing or decreasing to the programmed amplitude or until turned off.

Status row - The top row on the Therapy screen. This row displays if your therapy is on or off.

Stimulation - The delivery of electrical pulses to an appropriate area of the brain. Your therapy is actually mild electrical stimulation.

Stimulation cycling - This feature turns your stimulation on and off at regular intervals and is programmed by your doctor.

Synchronize - The process of sending and receiving information between the patient programmer and neurostimulator.

Therapy - Treatment of a disease or condition. When neurostimulation therapy is prescribed, a neurostimulation system is used to deliver stimulation to one or more areas of the brain.

Glossary

Therapy screen - The main screen displayed on the patient programmer.

Therapy settings - Your therapy can be adjusted by changing the rate, amplitude, or pulse width of the electrical stimulation. Your clinician programs all therapy settings. You may be able to adjust some therapy settings within clinician-defined limits. Available in Advanced mode only.

Warning - A statement describing an action or situation that could harm the patient.

Warning screen - A screen displayed on the patient programmer that alerts you to a problem with the programmer, antenna, or neurostimulator.

Glossary

16 English 37642 2020-01-01

M929532A050 Rev A 2020-01-01

Filename Date Time
UC200xxxxxx EN
4.6 x 6 inches (118 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.3: 09-13-2019

1 Introduction

M929532A059 Rev A 2020-01-01

How to use this manual

Refer to this manual after you receive an implanted neurostimulator. Ask your doctor to explain anything that is unclear.

- A glossary is provided at the beginning of this manual to describe terms that may be unfamiliar to you.
- Chapter 1 "Introduction", on page 17, describes how to use this manual, the purpose of the patient programmer, and provides general precautions related to using your patient programmer.
- Chapter 2 "Using your patient programmer", on page 25, describes the patient programmer and how to perform specific tasks.
- Chapter 3 "Adjusting therapy settings with Advanced mode", on page 65, explains

how to use your patient programmer to adjust your therapy settings.¹

- Chapter 4 "MRI examinations", on page 83, provides information about what you should do if you have an MRI examination.
- Chapter 5 "Important information about your rechargeable neurostimulator", on page 97 emphasizes the importance of checking and recharging your neurostimulator battery (for rechargeable devices only).
- Chapter 6 "Troubleshooting", on page 109, describes patient programmer warning and information screens, how to solve possible problems, and who to contact if your device is lost or broken.
- Chapter 7 "Maintenance", on page 127 describes how to care for your patient

¹ Your doctor has programmed your implanted neurostimulator to provide the correct therapy settings for your medical condition and will activate Advanced mode if appropriate.

programmer, including how to change the batteries.

For important safety information

Please see your DBS Patient Therapy Guide and DBS Therapy-specific Patient Booklet for important safety information, related warnings and precautions, and additional information about your therapy.

Purpose of your patient programmer

The Medtronic Model 37642 Patient Programmer is designed to operate with the following Medtronic neurostimulators:

- Activa™ PC Model 37601 (nonrechargeable battery)
- Activa RC Model 37612 (rechargeable battery)
- Activa SC Model 37602 (nonrechargeable battery)

Introduction 1

20 English 37642 2020-01-01

M929532A059 Rev A 2020-01-01

- Activa SC Model 37603 (nonrechargeable battery)

Purpose of the neurostimulation system (indications)

Refer to the indications sheet that is packaged with your patient programmer for the purpose of the neurostimulation system and related information.

Precautions

Patient control devices may affect other implanted devices—Do not place the patient control device (ie, patient programmer) over another type of active implanted medical device (eg, pacemaker, defibrillator, another type of neurostimulator). The patient control device could unintentionally change the operation of the other device.

Patient programmer handling—To avoid damaging the patient programmer, do not immerse the device in liquid; do not clean it

Introduction 1

37642 2020-01-01 English 21

with bleach, nail polish remover, mineral oil, or similar substances; and do not drop it or handle it in a way that might damage it.

Patient programmer use—The device is not certified for use in the presence of a flammable or anesthetic mixture with air or with oxygen or nitrous oxide. The consequences of using the device near flammable atmospheres are unknown.

Patient programmer modification—Do not modify this equipment. Modification of this equipment can result in damage to the programmer, causing the programmer to malfunction or become unusable.

Communication interference from EMI—When using your patient programmer to communicate with your neurostimulator, move away from equipment that may generate electromagnetic interference (EMI) or turn off the likely source of EMI. EMI may disrupt communication between the patient programmer and neurostimulator. Examples of EMI sources are computer monitors, cellular telephones, and motorized

Filename Date Time
UC200xxxxxx EN
4.6 x 6 inches (118 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.3: 09-13-2019

wheelchairs. For more information about EMI,
refer to your DBS Therapy-specific patient
booklet.

Introduction 1

37642 2020-01-01 English 23

M929532A059 Rev A 2020-01-01

Filename Date Time
UC200xxxxx EN
4.6 x 6 inches (118 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.3: 09-13-2019


Introduction 1

24 English 37642 2020-01-01

M929532A059 Rev A 2020-01-01

Filename Date Time
UC200xxxxxx EN
4.6 x 6 inches (118 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.3: 09-13-2019

2 Using your patient programmer

M929532A059 Rev A 2020-01-01

How the patient programmer works

The patient programmer is used to control and monitor your implanted neurostimulator or neurostimulators. For example, you will use your patient programmer to:

- Turn your therapy on or off.
- Check the neurostimulator and patient programmer battery status.
- Alert you when you need to check the status of the neurostimulator battery.¹
- Change therapy settings.²

Note: Make sure to keep your patient programmer with you at all times in the event that you need to adjust your therapy or turn your therapy off.

¹ Your clinician will turn this feature on or off and also set the time of day that it will alert you.

² Your clinician will set up the therapy settings you will be able to adjust depending on your specific therapy requirements. Discuss this with your doctor.

The patient programmer communicates with your neurostimulator by sending signals to and receiving signals from your implanted neurostimulator.

When you use the patient programmer, hold it directly over your implanted neurostimulator so that the programmer screen is facing out. The back of the patient programmer should be as close to the neurostimulator as possible, as shown in Figure 2.1.


Figure 2.1 Place the patient programmer directly over your neurostimulator.

A detachable antenna is available for use with your neurostimulator system. The detachable antenna makes it easier to view the patient programmer screen while adjusting your therapy settings or checking the neurostimulator battery status. See "Using the detachable antenna" on page 61 for more information.

Using your patient programmer 2

28 English 37642 2020-01-01

M029532A059 Rev A 2020-01-01

Multiple neurostimulators

Some patients have decided with their doctors to have two neurostimulators implanted.

If you have two neurostimulators implanted, you can use one patient programmer to communicate with both neurostimulators, but only with one neurostimulator at a time. The patient programmer only displays information about the neurostimulator with which it is currently synchronized.

The patient programmer must be turned off prior to synchronizing with the other neurostimulator or you may not be viewing the desired information.

Simple and advanced therapy modes

There are two therapy modes available for your patient programmer: Simple mode and Advanced mode. Your doctor will discuss which mode you will use.

- Simple mode uses the therapy settings selected by your doctor. Only your doctor can adjust your therapy settings.
- Advanced mode uses therapy settings selected by your doctor and also may allow you to change certain therapy settings and/or select preset therapy groups.

Using your patient programmer 2

30 English 37642 2020-01-01

M929532A059 Rev A 2020-01-01

Table 2.1 Tasks available for Simple and Advanced therapy modes

Tasks	Simple	Advanced
Check neurostimulator battery	Yes	Yes
Check patient programmer battery	Yes	Yes
Check status of neurostimulator (on or off)	Yes	Yes
Turn audio on or off. This does not change the Patient alert.	Yes	Yes
Adjust time	Yes	Yes
Adjust screen contrast	Yes	Yes
Adjust therapy settings	No	Yes ^a
Change active group	No	Yes ^a

^a This advanced feature is determined by your doctor.

For details about Advanced mode, see "Adjusting therapy settings with Advanced mode" on page 65.

Using your patient programmer 2

Summary of keys

Figure 2.2 Keys used to operate the patient programmer.

Refer to Table 2.2 for a description of the patient programmer keys.

Using your patient programmer 2

Table 2.2 Summary of keys

Key	Function
 Therapy On/Off	Turns your therapy on or off. - The patient programmer (or detachable antenna) must be held directly over the neurostimulator while pressing the Therapy On/Off key. - After using the Therapy On/Off key, verify your therapy status by viewing the displayed Therapy screen. Press this key to turn your therapy on or off.
 Check	Synchronizes the neurostimulator and patient programmer. Activates a selected group (Advanced mode only). Activates the "Return to Clinician Settings" (Advanced mode only). The patient programmer (or detachable antenna) must be held directly over the neurostimulator while pressing the Check  key. Press this key to synchronize the patient programmer and neurostimulator.

Using your patient programmer 2

Table 2.2 Summary of keys (continued)

Key	Function
 Power/ backlight On/Off	Turns the patient programmer and the backlight on or off. The backlight provides more light to the display. Press and release this key once to turn the patient programmer power on or off. Press and hold this key to turn the backlight on or off.

Table 2.2 Summary of keys (continued)

Key	Function
 Selection	Selects parameter setting for the left or right side of body (Therapy screen). Increases or decreases setting values (Parameter screen). When increasing or decreasing values, hold the patient programmer (or detachable antenna) directly over the neurostimulator and press the appropriate Selection key. You have two options to increase or decrease the settings with the Selection key. • Press and release the Selection key multiple times to increase or decrease the settings. • Press and hold the Selection key to increase or decrease the settings which will change approximately every half-second. Press the Selection key under an option displayed on the screen to select that option.

Using your patient programmer 2

Table 2.2 Summary of keys (continued)

Key	Function
	Moves the selection box up or down on the Therapy screen or to the next available screen.
Navigator	The Options  icon at the end of a row on the Therapy screen indicates that the row continues on the next screen. Press the left  or right  arrows on this key to navigate to the next screen. Press the up  or down  arrows on this key to move the selection box to the next row.

Using the Navigator key

The **Navigator** key moves the selection box on the patient programmer screen to the right or left and up or down (Figure 2.3).

Figure 2.3 Navigator key.

- To move the selection box across a row that continues onto the next screen, press the left ◀ or right ▶ arrows on the **Navigator** key.
- To move the selection box to a new row, press the up ▲ and down ▼ arrows on the **Navigator** key.

The image below is an example of how the **Options** icon is displayed. The **Options** icon is displayed at the end of a row when that row has additional items to view to the right or left (Figure 2.4).


Figure 2.4 The **Options** icon and selection box.

When moving the selection box with the **Navigator** key, you do not need to hold the patient programmer over your neurostimulator.

Using your patient programmer 2

The Therapy screen

The **Therapy** screen displays your neurostimulator settings and battery status.

- In Simple mode:

The **Therapy** screen shows whether your neurostimulator is on or off and the status of the neurostimulator battery. The Simple mode screen appears the same for Activa PC, RC, and SC models (Figure 2.5).


Figure 2.5 The **Therapy** screen in Simple mode for Activa PC, RC, and SC models

- In Advanced mode:

The **Therapy** screen shows whether your neurostimulator is on or off, the battery status, and the specific parameters set for your therapy.

In most cases, the Advanced mode screen appears differently for Activa SC models than the Activa PC and Activa RC models (Figure 2.6). In most cases, one Activa PC or RC neurostimulator stimulates two sides of the brain. In contrast, one Activa SC neurostimulator stimulates one side of the brain.


Figure 2.6 The Therapy screen in Advanced modes Activa PC, RC, and SC models.

Using your patient programmer 2

Table 2.3 Therapy screen icons

Row	Icons	Description
Status		Therapy is on
		Therapy is off
Battery (Activa RC)		Neurostimulator battery charge level (rechargeable) ^a
	 Needs Charge Full Charge	
Battery (Activa PC and SC models)		Neurostimulator battery status (nonrechargeable)
	 OK  ERI Battery OK Battery near end of service	
Parameter/ group ^b		Active group ^c
		Amplitude
		Pulse width
		Rate

^a For complete rechargeable battery information, see Chapter 5 "Important information about your rechargeable neurostimulator".

^b The **Parameter/group** row only appears in Advanced mode.

Using your patient programmer 2

Table 2.3 Therapy screen icons

Row	Icons	Description
-----	-------	-------------

^c If only one group is available, this row displays only the parameter setting icon.

Synchronizing the patient programmer and your neurostimulator

In order to check your neurostimulator status, change therapy settings, or turn your neurostimulator on and off, you must first synchronize the patient programmer and the neurostimulator.

Complete the following steps to synchronize the patient programmer and neurostimulator.

1. Locate the **Check** ☒ key.
2. Hold the patient programmer (or antenna) directly over the neurostimulator with the screen facing outward.
3. Press the **Check** ☒ key.

Using your patient programmer 2

Note: Pressing the **Check**  key also turns on the patient programmer (Figure 2.7).

Use the
Check  key
to synchronize

Figure 2.7 Press the **Check** key to synchronize your neurostimulator and patient programmer.

4. After pressing the **Check**  key, the **Communication** screen appears briefly.


Figure 2.8 Communication screen.

Using your patient programmer 2

- a. If synchronization is successful, the Therapy screen appears (Figure 2.9).

Figure 2.9 Therapy screen (Simple mode).

- b. If synchronization is not successful, the Position antenna screen or the Poor communication screen appears (Figure 2.10 and Figure 2.11). See Table 6.3 on page 115 for information on how to proceed when these screens appear.


Figure 2.10 Position antenna screen


Figure 2.11 Poor communication screen.

Using your patient programmer 2

Note: See Chapter 6 "Troubleshooting" for more information about synchronizing the patient programmer and neurostimulator.

Turning your therapy on or off

Complete the following steps to turn your therapy on. Be sure to talk to your doctor about turning your therapy on or off before following this procedure.

1. Synchronize the patient programmer and neurostimulator.
 - a. Hold the patient programmer (or antenna) directly over your neurostimulator with the screen facing outward.
 - b. Press the Check  key. The Therapy screen appears.

The **Status** icon on the Therapy screen indicates whether your therapy is turned on or off (Figure 2.12).

Using your patient programmer 2

Figure 2.12 Therapy status is displayed.

2. Press the **Therapy On/Off** key (⚡).

If your therapy is currently off, pressing the **Therapy On/Off** key will turn the therapy on.

If your therapy is currently on, pressing the **Therapy On/Off** key will turn the therapy off.

The display flashes **Off** when the therapy is turned off. If the display is set to icon mode, a warning symbol (⚠) flashes when the therapy is turned off.


Figure 2.13 Therapy off in icon mode.

Using your patient programmer 2

Note: See "Changing preferences" on page 51 for more information about icon mode.

The neurostimulator battery

Your neurostimulator battery may be either a rechargeable battery or a nonrechargeable battery. Check with your doctor if you do not know what type of neurostimulator battery you have.

About rechargeable batteries

The rechargeable neurostimulator battery needs to be checked daily to make sure the charge level is adequate for your therapy.

Note: For complete information about the rechargeable neurostimulator battery, see Chapter 5 "Important information about your rechargeable neurostimulator" on page 97.

Using your patient programmer 2

46 English 37642 2020-01-01

M929532A059 Rev A 2020-01-01

Checking the nonrechargeable battery status

1. Synchronize the patient programmer and neurostimulator.

Note: For information on checking and charging the rechargeable neurostimulator battery, see Chapter 5 "Important information about your rechargeable neurostimulator".

- a. Hold the patient programmer (or antenna) directly over your neurostimulator with the screen facing outward.
- b. Press the **Check**  key. The **Therapy** screen appears (Figure 2.14).

 On	 Off
Therapy is turned On	Therapy is turned Off

Figure 2.14 Therapy status is displayed.

Using your patient programmer 2

2. Check the current status reading for the battery:

- Simple Mode: Press the left ◀ or right ▶ arrow on the **Navigator** key to display the **Battery status** screen (Figure 2.15).
- Advanced mode: press the up ▲ or down ▼ arrow on the **Navigator** key to select the **Battery** row. Press the left ◀ or right ▶ arrow on the **Navigator** key to display the **Battery status** screen (Figure 2.15).

Figure 2.15 Battery Status screen.


Possible nonrechargeable battery status conditions

The patient programmer displays an information screen if the neurostimulator battery needs your attention (Table 2.4).

Using your patient programmer 2

See Chapter 6 "Troubleshooting" for more information about battery-related screens.

Table 2.4 Nonrechargeable neurostimulator battery screens

	If the Therapy screen displays ERI, the neurostimulator battery is close to its end of service date.^a You should call your doctor to report.
	This information screen appears when you synchronize your neurostimulator and patient programmer when the nonrechargeable battery status is low. This means that your therapy will not be available soon. Press any key to clear this message. Call your doctor to report this message.

Using your patient programmer 2

Table 2.4 Nonrechargeable neurostimulator battery screens (continued)

EOS condition

If the EOS screen displays, then the battery is at end of service (EOS), therapy has stopped, and the neurostimulator needs to be replaced.

Call your doctor immediately.

^a If your nonrechargeable neurostimulator battery is approaching ERI condition, your doctor may set up an alert on the patient programmer to remind you to check your neurostimulator battery voltage once a day. For more information see, "Patient programmer alert" on page 57.

Using your patient programmer 2

Customizing your patient programmer

You can adjust the audio and display settings for your patient programmer by changing the preferences. Preferences can be reached from the **Status** row of the **Therapy** screen. Changing the preferences for the patient programmer do not affect your therapy.

See Table 2.5 for a list of the available preference settings.

Table 2.5 Preference icons

Icons	Preference
	Audio
	Contrast
	Time ^a
	Time and number format
Abc...  	Text or icon only display

^a The Patient alert feature is based on the time set in your patient programmer. For more information, see "Patient programmer alert" on page 57.

Changing preferences

Complete the following steps to adjust the patient programmer preferences.

1. Synchronize the patient programmer and neurostimulator.
 - a. Hold the patient programmer (or antenna) directly over your

Using your patient programmer 2

neurostimulator with the screen facing outward.

- b. Press the Check  key. The Therapy screen appears (Figure 2.16).

Figure 2.16 Access preferences from the Status row.

2. If in Advanced mode, press the down  arrow on the Navigator key to move the selection box to the Status row.
3. Press the left  or right  arrow on the Navigator key to move the selection box to the desired preference (Figure 2.17).

Using your patient programmer 2

Figure 2.17 Desired preference (audio in this example).

4. Press the down ▼ arrow on the **Navigator** key to move the selection box to the **Change** row (Figure 2.18).



Figure 2.18 Change row (audio in this example).

5. Follow the steps in Table 2.6 to change the selected preference.

Table 2.6 Changing programmer preferences

Audio preference

- Press the left ◀ or right ▶ arrow on the Navigator key to move the selection box to audio on ⏮ or off ⏭.
- Continue to step 6 through step 8 on page 56.

Note: The Patient alert will play even if the Audio preference is turned off.


Contrast preference

- Press the left ◀ or right ▶ arrow on the Navigator key to move the selection box to the contrast slider bar.
- Press the Selection keys to make the contrast lighter ○ or darker ●.
- Continue to step 6 through step 8 on page 56.


Using your patient programmer 2

**Table 2.6 Changing programmer preferences
(continued)**

Time preference

- Press the left ◀ or right ▶ arrow on the **Navigator** key to move the selection box to the hour, minutes, or time of day (A or P).
- Press the **Selection** keys to increase or decrease the displayed time.
- Continue to step 6 through step 8 on page 56.



Time and number format preference

- Press the left ◀ or right ▶ arrow on the **Navigator** key to move the selection box to a 12-hour clock or to a 24-hour clock.
- Continue to step 6 through step 8 on page 56.

Using your patient programmer 2

**Table 2.6 Changing programmer preferences
(continued)**

Text on	Text or icon display format
	a. Press the left ◀ or right ▶ arrow on the Navigator key to move the selection box to the Text on or Text off/icon only display setting.
Text off/icon only 	b. Continue to step 6 through step 8 on page 56. Note: The text displayed on the patient programmer will always be in English.

Using your patient programmer 2

- Press the up ▲ arrow on the **Navigator** key to return the selection box to the **Status** row.
- Press the left ◀ or right ▶ arrow on the **Navigator** key to return to the **Therapy** screen.
- Hold the patient programmer or antenna directly over your neurostimulator and press the **Check** ✓ key to save the preference change.

Patient programmer alert

Your doctor may set up an alert on the patient programmer to remind you to check your neurostimulator battery once a day (rechargeable neurostimulator batteries only).

Note: For nonrechargeable neurostimulator batteries, your doctor may set up this alert for you when your battery is near its end of service.

The default setting for the alarm is 11:00 a.m. However, your doctor may change this to a different time of day if appropriate.

The alert time is based on the time set for the patient programmer. To adjust this time, see Table 2.6 on page 54.

The patient programmer alert beeps for 30 seconds, with a pause every 10 seconds. The **Patient alert** screen displays while the alarm is beeping (see Figure 2.19). The Patient alert sounds even if the audio preference is turned off.

Using your patient programmer 2


Figure 2.19 Patient alert screen.

This alert repeats every 15 minutes until the neurostimulator battery status has been checked.

To check battery status, see "Checking the nonrechargeable battery status" on page 47 or "Checking the rechargeable neurostimulator battery" on page 102.

Patient programmer accessories

Carrying case and patient programmer ID label

A carrying case and identification label are included with your patient programmer. Two AAA alkaline batteries that provide the power for your patient programmer are also included.

The carrying case has a pouch to hold the patient programmer and the quick reference guide.

Figure 2.20 Insert the patient programmer into the case.

You can attach the patient programmer to a belt using the loop on the back of the case.

Using your patient programmer 2

Place the provided identification label on the back of your patient programmer just in case the patient programmer is lost (Figure 2.21).


Figure 2.21 Place the adhesive label on the back of the patient programmer.

Using your patient programmer 2

Detachable antenna

A detachable antenna (Model 37092) is useful for viewing the patient programmer screen while you are adjusting therapy. When the antenna is connected to the patient programmer, hold the antenna or attach it to your clothing directly over the neurostimulator to adjust settings or check the neurostimulator battery.

Using the detachable antenna

1. Place the antenna directly over your neurostimulator (Figure 2.22).


Figure 2.22 Place the antenna directly over your neurostimulator.

2. Pull the fabric of your clothing through the large opening in the antenna. Then, wedge the fabric in the narrow slit to secure the antenna in place (Figure 2.23).

Using your patient programmer 2


Figure 2.23 Pull the fabric through the slit (a) and wedge in place (b).

3. Push the antenna plug firmly into the antenna jack (Y) on the patient programmer (Figure 2.24).

Figure 2.24 Insert the antenna plug into the antenna jack.


Using your patient programmer 2

After the antenna is connected, follow the instructions for using the patient programmer.

When you have finished using the patient programmer, grasp the antenna plug and pull it out.

△ **Caution:** Do not pull directly on the antenna cable to disconnect the cable from the programmer because this may damage the antenna cable.


Using your patient programmer 2

Filename Date Time
UC200xxxxxx EN
4.6 x 6 inches (118 mm x 152 mm)


Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.3: 09-13-2019


Using your patient programmer 2

64 English 37642 2020-01-01

M929532A050 Rev A 2020-01-01

Filename Date Time
UC200xxxxx EN
4.6 x 6 inches (118 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.3: 09-13-2019

3 Adjusting therapy settings with Advanced mode

M929532A059 Rev A 2020-01-01

Adjusting your therapy settings

You may find that on some days, in order to receive the most effective therapy, you need to adjust your therapy several times. On other days you may not need to adjust it at all. Your doctor will provide complete guidelines about when you may want to adjust your therapy.

- Ask your doctor to print a report with your programmed settings.
- When a therapy setting is changed, you will see the change on the **Therapy** screen.
- If your patient programmer audio is turned on, you will hear one tone that means the change was effective.
- Three rapid tones means that the change could not be confirmed or that the synchronization was not successful.
- Verify the current therapy settings by pressing the **Check**  key and viewing the **Therapy** screen.

Tips for adjusting your therapy:

- Always adjust the therapy settings to the lowest possible setting that provides symptom relief.
 - High therapy settings use more rechargeable battery charge which may force you to charge the battery more frequently.
 - High therapy settings may shorten the service life of the nonrechargeable battery.
- Before changing any therapy settings with your patient programmer, you need to synchronize the patient programmer with the neurostimulator.

Note: If you have two neurostimulators implanted, you can use one patient programmer to communicate with both neurostimulators, but only with one neurostimulator at a time. To communicate with another neurostimulator, turn the patient programmer off and then press the **Check**

Adjusting therapy settings with Advanced mode 3

 key to synchronize while holding over the desired neurostimulator.

- If you experience uncomfortable therapy or unwanted side effects when increasing the amplitude, pulse width, or rate, decrease the setting to a comfortable therapy setting. You may also use the **Therapy On/Off**  key to stop the therapy.

Adjusting therapy settings with Advanced mode 3

Therapy settings

If your neurostimulator is set to Simple mode, you will not be able to adjust your therapy settings. Your doctor can adjust your settings when necessary.

Notes:

- To increase amplitude, pulse width, or rate, your therapy must be turned on.
- To decrease amplitude or pulse width, your therapy may be turned on or off.
- To decrease rate, your therapy must be turned on.

- For each group, only one parameter can be changed. For example, group A may allow you to change only the amplitude of the therapy and group B may allow you to change only the rate of the therapy.
- A group may also be programmed for view only and you will not be able to adjust any settings. Discuss your settings with your doctor. For more information about groups, see "Selecting a new group" on page 73.

Changing therapy settings in Advanced mode

The ability to change therapy settings is only available in Advanced mode.

You can only change the settings of the active group. If you want to change settings available for a non-active group, you will need to make that group active first. See "Selecting a new group" on page 73 for complete instructions.

Adjusting therapy settings with Advanced mode 3

1. Synchronize the patient programmer and neurostimulator.
 - a. Hold the patient programmer (or antenna) directly over your neurostimulator with the screen facing outward.
 - b. Press the **Check** (✓) key. The **Therapy** screen appears.
2. Press the down ▼ arrow on the **Navigator** key to move the selection box to the **Parameter** row (Figure 3.1).

Adjusting therapy settings with Advanced mode 3



Figure 3.1 Parameter row for Activa PC, RC, and SC models.

3. Press the **Selection** key directly under the displayed setting you want to adjust (Figure 3.2).

Parameter row for Activa PC and RC Two sides of the brain

Parameter row for Activa SC Models One side of the brain


Figure 3.2 Select the appropriate Selection key
on Activa PC, RC, and SC models.

Adjusting therapy settings with Advanced mode 3

4. Press the **Selection** keys again to increase or decrease the selected therapy setting as needed (Figure 3.3).¹


Figure 3.3 Increasing or decreasing a therapy setting.

Notes:

- Hold the patient programmer or detachable antenna over your neurostimulator when increasing or decreasing your settings. The therapy settings will immediately change as you press the **Selection** keys.

¹ If you attempt to increase or decrease the therapy setting beyond the available limits, an information screen appears, see Table 3.1 on page 81.

- Press and hold down the **Selection** key to increase or decrease the setting which will change approximately every half-second.

5. Press any arrow on the **Navigator** key to return to the **Therapy** screen.
6. Repeat steps 2-5 for the other side of your body if desired.

Note: If you have multiple neurostimulators, the patient programmer must be turned off prior to synchronizing with the other neurostimulator.

Selecting a new group

The group feature is only available in Advanced mode.

△ **Caution:** Select the group that your clinician has recommended for your specific needs.

1. Synchronize the patient programmer and neurostimulator.

Adjusting therapy settings with Advanced mode 3

- a. Hold the patient programmer (or antenna) directly over your neurostimulator with the screen facing outward.
 - b. Press the Check ☒ key. The Therapy screen appears.
2. Press the down ▼ arrow on the **Navigator** key to move the selection box to the **Parameter** row (Figure 3.4).

Adjusting therapy settings with Advanced mode 3


Figure 3.4 Parameter row.

3. Press the right ► arrow on the **Navigator** key to display the available groups. (Figure 3.5).


Figure 3.5 Group screen.

4. Press the up ▲ or down ▼ arrow on the Navigator key to move the selection box the desired group (Figure 3.6).

Figure 3.6 Select the desired group.

5. Resynchronize the patient programmer and neurostimulator.
 - a. Hold the patient programmer (or antenna) directly over your neurostimulator with the screen facing outward.
 - b. Press the Check  key to synchronize with the patient programmer and activate the selected group.

Adjusting therapy settings with Advanced mode 3

6. The **Communication** screen briefly appears, followed by the **Therapy** screen (Figure 3.7).


Figure 3.7 Communication screen.

The **Therapy** screen will display your new therapy parameters and group setting (Figure 3.8).

Adjusting therapy settings with Advanced mode 3

Figure 3.8 Therapy screen with new group setting.

Returning your therapy to original settings

If you have changed your therapy settings and want to return to the settings selected by your doctor, follow these steps:

1. Synchronize the patient programmer and neurostimulator.
 - a. Hold the patient programmer (or antenna) directly over your neurostimulator with the screen facing outward.
 - b. Press the **Check**  key. The **Therapy** screen appears.
2. Press the down  arrow on the **Navigator** key to move the selection box to the **Parameter** row (Figure 3.9).

Adjusting therapy settings with Advanced mode 3

Figure 3.9 Parameter row.

3. Press the right ► arrow on the **Navigator** key to display the available groups. The **Group settings** screen appears (Figure 3.10).


Figure 3.10 Group settings screen.

To change from the current group settings back to the doctor default settings:

- a. With the selection box highlighting the current group, press the right ► arrow on the **Navigator** key. The **Group reset** screen appears (Figure 3.11).


Figure 3.11 Group reset screen.

- b. Press the Check ☒ key.

To change an inactive group back to the doctor default settings:

- a. Press the up  or down  arrow on the Navigator key to move the selection box to the desired inactive group (Figure 3.12).


Figure 3.12 Select new group.

- b. Press the right  arrow on the Navigator key. The Group options screen appears (Figure 3.13).

Adjusting therapy settings with Advanced mode 3


Figure 3.13 Group options screen.

- c. Press the down ▼ arrow on the **Navigator** key to move the selection box to the **Clinician setting** option (Figure 3.14).


Figure 3.14 Change group to clinician settings.

- d. Press the right ► arrow on the **Navigator** key to display the **Group reset** screen (Figure 3.15).


Figure 3.15 Group reset screen

- e. Press the **Check** ☒ key.



Adjusting therapy settings with Advanced mode 3

4. The patient programmer beeps when the group settings are changed (if the audio is turned on).

Understanding therapy parameter limits


Each therapy parameter will have a default limit set by your doctor. If you attempt to change a therapy parameter outside these set limits, you will see one of the screens shown in Table 3.1.

Table 3.1 Parameter limit screens

Lower limit	
	You tried to decrease a parameter (amplitude, pulse width, or rate) below the lowest value allowed. Press any arrow on the Navigator key to clear the screen.
	
	


Adjusting therapy settings with Advanced mode 3

**Table 3.1 Parameter limit screens
(continued)**

Upper limit	You tried to increase a parameter (amplitude, pulse width, or rate) above the highest value allowed. Press any arrow on the Navigator key to clear the screen.
	
	
	

Adjusting therapy settings with Advanced mode 3

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 5.3: 09-13-2019

4 MRI examinations

 ООО ЦПМИ  INFO@CIRMI.RU  WWW.GOSZDRAVNADZOR.RU

If you have an MRI appointment

 **MR Conditional**—Depending on what kind of neurostimulation system components that you have implanted, you may be eligible for one of the following types of magnetic resonance imaging (MRI) scans:

- MRI scans of the head only (this is referred to as head-only eligible).
- MRI scans of any part of your body (this is referred to as full-body eligible).

Note: You should inform your clinician managing your neurostimulation system that an MRI examination has been prescribed for you and that you need him or her to determine what type of MRI scan you are eligible to receive.

MRI examinations 4

Responsibilities of the patient in preparing for the MRI appointment

Bring the following to every MRI appointment:

- Your patient identification (ID) card for your neurostimulation system if you received one from Medtronic.

Note: If you have two neurostimulation systems implanted in your body, bring both ID cards to your MRI appointment.

- Your MRI scan eligibility sheet if you were given one from your clinician managing your neurostimulation system.
- Your patient programmer.

For more information about how MRI can affect your neurostimulation system and what you should do if you have an MRI appointment, refer to the *DBS Patient Therapy Guide*.

Preparing your neurostimulation system for the MRI scan

Your implanted neurostimulation system (ie, therapy) may need to be turned off prior to your MRI scan or reprogrammed to allow the

MRI examinations 4

system to safely remain on during your MRI scan. This will depend on the neurostimulator model implanted in your body, the therapy settings of your neurostimulator, and the type of MRI scan you are eligible to receive.

Your clinician managing your neurostimulation system will inform you whether your system should be on or off during the MRI scan.

If therapy should be turned off prior to your MRI scan, refer to "Turning off therapy before the MRI scan" on page 86.

If therapy can remain on during your MRI scan, refer to "Activating settings to allow therapy on during the MRI scan" on page 90.

MRI examinations 4

Turning off therapy before the MRI scan

If your clinician managing your neurostimulation system indicates that therapy should be off during the MRI scan, therapy can be turned off with the patient

programmer (see page 87) or the clinician programmer (see page 88).

Using the patient programmer to turn off therapy

If you brought your patient programmer to the MRI appointment, you can turn off therapy before your MRI scan and outside of the MRI scanner (magnet) room.

For instructions on turning off therapy, go to "Turning your therapy on or off" on page 44.

During the MRI scan, keep therapy off.

⚠ Caution: Do not turn therapy back on before your MRI scan. Leaving therapy on during the scan could increase the potential for unintended stimulation.

After the MRI scan and outside of the MRI scanner (magnet) room, you can turn therapy back on using your patient programmer or you can return to the clinician managing your neurostimulation system to have therapy turned back on (refer to "Turning therapy back on after the MRI scan" on page 89).

MRI examinations 4

Note: Do not take the patient programmer into the MRI scanner (magnet) room. If you brought the recharger with you to the MRI appointment, do not take the recharger into the MRI scanner (magnet) room either.

Using the clinician programmer to turn off therapy

Your clinician managing your neurostimulation system can turn off therapy before your MRI examination using the clinician programmer.

Do not turn therapy back on before your MRI scan.

Caution: Do not turn therapy back on before your MRI scan. Leaving therapy on during the scan could increase the potential for unintended stimulation.

The clinician may give you an MRI scan eligibility sheet to bring to your MRI appointment. Give the eligibility sheet to the MRI clinician.

During the MRI scan, keep therapy off.

MRI examinations 4

After the MRI scan and outside of the MRI scanner (magnet) room, you can turn therapy back on using your patient programmer or you can return to the clinician managing your neurostimulation system to have therapy turned back on (refer to "Turning therapy back on after the MRI scan" on page 89).

Turning therapy back on after the MRI scan

Turn therapy back on when the MRI scan is complete and you are **outside** of the MRI scanner (magnet) room.

Therapy can be turned on with the patient programmer or the clinician programmer:

- If you do not have your patient programmer with you, go to the clinician managing your neurostimulation system to turn on therapy using the clinician programmer.
- If you brought your patient programmer to the MRI appointment, turn therapy back on using the patient programmer. For

MRI examinations 4

instructions on turning on therapy, go to
"Turning your therapy on or off" on
page 44.

Note: If a power-on-reset (POR) screen
appears on the patient programmer, see
Table 4.1 on page 95.

Activating settings to allow therapy on during the MRI scan

Depending on your therapy settings, your
clinician managing your neurostimulation
system may create a new group that should
be activated prior to your MRI scan. This will
allow your system to safely remain on during
the scan.

If your clinician managing your
neurostimulation system creates a new group
setting so that therapy can remain on during
your MRI scan, the new group can be
activated with the patient programmer (see
page 91) or the clinician programmer (see
page 92).

Using the patient programmer to activate a new group

If you brought your patient programmer to the MRI appointment, you can activate the new group before your MRI scan and outside of the MRI scanner (magnet) room. This will allow you to keep therapy on during your MRI scan.

For instructions on activating a new group, go to "Selecting a new group" on page 73.

During the MRI scan, keep therapy on, but do not return therapy to your original group setting.

⚠ **Caution:** Do not return therapy to your original group setting before your MRI scan. Returning therapy to your original group setting could increase the potential for unintended stimulation during the scan.

After the MRI scan and outside of the MRI scanner (magnet) room, you can return therapy to your original group setting using

MRI examinations 4

your patient programmer or you can return to the clinician managing your neurostimulation system to reprogram therapy to your original group setting (refer to "Returning therapy to your original group setting after the MRI scan" on page 94).

Note: Do not take the patient programmer into the MRI scanner (magnet) room. If you brought the recharger with you to the MRI appointment, do not take the recharger into the MRI scanner (magnet) room either.

Using the clinician programmer to activate a new group

The clinician managing your neurostimulation system can activate the new group before your MRI examination using the clinician programmer. This will allow you to keep therapy on during your MRI scan.

Do not change therapy to your original group setting before the MRI scan.

△ **Caution:** Do not return therapy to your original group setting before your MRI

MRI examinations 4

scan. Returning therapy to your original group setting could increase the potential for unintended stimulation during the scan.

The clinician may give you an MRI scan eligibility sheet to bring to your MRI appointment. Give the eligibility sheet to the MRI clinician.

During the MRI scan, keep therapy on, but do not return therapy to your original group setting.

After the MRI scan and outside of the MRI scanner (magnet) room, you can return therapy to your original group setting using your patient programmer or you can return to the clinician managing your neurostimulation system to reprogram therapy to your original group setting (refer to "Returning therapy to your original group setting after the MRI scan" on page 94).

Returning therapy to your original group setting after the MRI scan

Return therapy to your original group setting when the MRI scan is complete and you are **outside** of the MRI scanner (magnet) room.

Group settings can be changed with the patient programmer or the clinician programmer:

- If you do not have your patient programmer with you, go to the clinician managing your neurostimulation system to return therapy to your original group setting using the clinician programmer.
- If you brought your patient programmer to the MRI appointment, return therapy to your original group setting using the patient programmer. For instructions on returning therapy to your original group setting, go to "Selecting a new group" on page 73.

MRI examinations 4

94 English 37642 2020-01-01

M929532A059 Rev A 2020-01-01

Note: If a power-on-reset (POR) screen appears on the patient programmer, see Table 4.1.

Table 4.1 POR screens




Screen	Description and action
	Error code POR: Your therapy has stopped. Call your clinician to restart your therapy.
Warning POR (hazard triangle symbol in the upper left corner)	

Table 4.1 POR screens (continued)

Screen	Description and action
 Informational POR (small "i" icon in the upper left corner)	The neurostimulator has been reset. Your therapy has stopped. 1. Press any arrow on the Navigator key to clear the screen 2. Reset the time, if prompted (for instructions, see Table 2.6). Before resetting the time, confirm the patient programmer or antenna is directly over the neurostimulator. 3. Press the up arrow on the patient programmer. This action will clear the POR. 4. Press the Therapy on/off  key to turn your therapy on. Call your clinician to report the reset.

Filename Date Time
UC200xxxxxx EN
4.6 x 6 inches (118 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.3: 09-13-2019

**5 Important
information about
your rechargeable
neurostimulator**

M929532A059 Rev A 2020-01-01

The rechargeable neurostimulator battery

The following chapter explains the importance of maintaining the battery charge level of your rechargeable neurostimulator.¹

⚠ **Warning:** It is very important to check your neurostimulator battery every day. If the therapy provided by your neurostimulator should stop abruptly this could cause your symptoms to return. In some cases, your symptoms may return at a greater intensity than before your implant. In rare situations, this could result in a medical emergency.

It is important for you to recharge your battery on a regular, frequent basis as recommended by your doctor to make sure that the neurostimulator battery remains charged. Depending on your therapy settings this may mean charging your battery once a day or once a week,

¹ This information only applies to implantable neurostimulators that use a rechargeable battery.

Important information about your rechargeable neurostimulator 5

37642 2020-01-01 English 99

Check and charge: make it a habit

Because this is so important, make it a priority to check and charge your neurostimulator battery on a regular schedule. See Table 5.1 on page 107 for rechargeable battery alert screens that indicate a low battery charge level.

Recommendations:

- Check your battery charge level at the same time every day (your doctor can set the patient programmer alert for this time).
- Combine checking your battery charge level with something else you do every day in order to make it a convenient habit.
- Make sure to bring your neurostimulator recharging system with you when you travel or are hospitalized (even for overnight).
- Allow enough time to fully charge the neurostimulator. Depending on the charge level of the battery when you begin

Note: For more information, see the manual packaged with your recharging device.

- **Make sure that you use your charging system as directed in your recharger manual. If the devices are not aligned correctly or if there is open space between the devices, it may take longer to charge the neurostimulator battery.**
- **Charge your neurostimulator battery before it is depleted. The battery should be charged before it has only 25% power remaining.**

You do not need to wait until the battery charge level is low. If it's more convenient, you can charge the battery every day.

Note: Remember that even when your therapy is turned off, the neurostimulator battery charge will still deplete over time. You

Important information about your rechargeable neurostimulator 5

should continue to check the battery daily and charge it when necessary.

Consult with your doctor about how often you should charge your neurostimulator battery based on your individual therapy settings.

Checking the rechargeable neurostimulator battery

1. Synchronize the patient programmer and neurostimulator.
 - a. Hold the patient programmer (or antenna) directly over your neurostimulator with the screen facing outward.

Press the Check ☒ key. The Therapy screen appears displaying the neurostimulator battery charge level (Figure 5.1).


Figure 5.1 The neurostimulator battery charge level is displayed on the Therapy screen.

Low battery status

△ Caution: Charge the neurostimulator when you see a **Low battery** screen on the patient programmer. If **LOW** is flashing, the battery charge level is at 0 and your therapy will stop soon unless you recharge the battery.

If you see the following screen (Figure 5.2), the neurostimulator battery is in a discharged state and your therapy has stopped.

Important information about your rechargeable neurostimulator 5


Figure 5.2 Neurostimulator battery is in a discharged state.

You should charge the battery as soon as possible. If the battery remains in a discharged state, it continues to lose charge and it will become overdischarged.

If the neurostimulator battery is allowed to overdischarge, you may not be able to recharge or communicate with the neurostimulator; however, your doctor may be able to restore the battery function.

Allowing the neurostimulator battery to overdischarge will permanently affect the neurostimulator in one of the following ways:

- Battery function is restored; however, charging sessions may be more

frequent because battery capacity has been reduced.

- Battery function is not restored and the neurostimulator must be surgically replaced. Battery function is not restored because:
 - the neurostimulator battery is permanently damaged.
 - the neurostimulator battery has been overdischarged and restored twice before. The third time the battery is overdischarged, the neurostimulator will reach end of service. Surgery is required to replace the neurostimulator.

Rechargeable battery overdischarge

You can avoid the consequences of an overdischarged battery and the loss of therapy by keeping the battery charged.

Important information about your rechargeable neurostimulator 5

Charging the neurostimulator battery




When the neurostimulator battery charge level is low, charge the battery as described in the manual packaged with the charging system.

Note: You can charge your battery more often, even daily. You do not need to wait until the battery charge level is low.

Your neurostimulator battery can be charged many times; however, eventually the neurostimulator will need to be replaced.

If the neurostimulator battery needs charging, you will see one of the screens shown in Table 5.1.

Table 5.1 Rechargeable neurostimulator battery screens

	The neurostimulator battery charge level is low. When the "Low" message is flashing, therapy will not be available soon.
 Charge battery screen	Charge your neurostimulator battery. Refer to the manual packaged with the charging system. Press any arrow on the Navigator key to clear this message from the screen.
 Battery is in a discharged state	The neurostimulator battery charge level is depleted and needs immediate charging. The battery is in a discharged state and therapy has stopped when the battery is depleted. Charge the neurostimulator battery now. Refer to the manual packaged with the charging system.^a

Important information about your rechargeable neurostimulator 5

Table 5.1 Rechargeable neurostimulator battery screens (continued)

EOS condition

If this screen displays, the battery is at the end of service, therapy has stopped, and the neurostimulator needs to be replaced.

If the EOS screen displays, call your doctor immediately.

^a If necessary, you can stop charging the battery before it is completely recharged in order to turn the therapy back on. Finish charging the battery after your therapy is turned back on.


Patient programmer alert

Your doctor may set up an alert on the patient programmer to remind you to check your neurostimulator battery once a day.

For more information, see "Patient programmer alert" on page 57.

Filename Date Time
UC200xxxxxx EN
4.6 x 6 inches (118 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.3: 09-13-2019

6 Troubleshooting

M929532A059 Rev A 2020-01-01

This chapter will help you solve problems with your patient programmer. It also provides information on when to call your doctor.

Note: If you cannot solve a problem or if your problem is not described here, contact your doctor.

Patient Programmer screens

After you synchronize your patient programmer with your implanted neurostimulator, the patient programmer may display a Warning , Communication , or Information  screen to provide you with information about your system, alert you to a problem with your system, or to guide you while using the patient programmer.



Warning screens

Warning screens indicate a problem with the patient programmer, the antenna, or the neurostimulator. If the patient programmer audio preference is turned on, three tones

alert you when a warning screen is displayed on the patient programmer.




Table 6.1 describes warning screens and provides instructions (**see blue text**) on how to resolve the problem and clear the screen.

Table 6.1 Warning screens

Screen	Cause and action
 Replace patient programmer batteries	The patient programmer batteries are depleted. Programming is not possible. Replace the patient programmer batteries now.
 Device not supported	The implanted device that you are attempting to communicate with is not compatible with the patient programmer. Call your doctor.
 Synchronize patient programmer and neurostimulator	The patient programmer and the neurostimulator must be synchronized. Synchronize the patient programmer and neurostimulator.


Troubleshooting 6

Table 6.1 Warning screens (continued)

Screen	Cause and action
Call your doctor screens (compare to the information screens on page 114.)	The system is not working correctly. Therapy may have stopped. Read the error code at the bottom of the screen.
	Error codes 0 to 252: Remove batteries from the patient programmer, wait several seconds, then re-insert the batteries. If the error message appears again, call your doctor.
	
	Error code EOS: Your neurostimulator has reached its end of service. Therapy is not available. Call your doctor.
	Error code POR: Your therapy has stopped. Call your doctor to restart your therapy.
	Other codes: Write down the code shown on the screen. Call your doctor.

Troubleshooting 6

Table 6.1 Warning screens (continued)

Screen	Cause and action
Rechargeable battery only  Rechargeable battery is in a discharged state	The rechargeable neurostimulator battery is in a discharged state and therapy has stopped. Charge the neurostimulator battery now. Refer to the manual packaged with the charging system.

Warning: If the neurostimulator battery is allowed to overdischarge, you may not be able to recharge it. However, your doctor may be able to restore the battery function.

Refer to "The rechargeable neurostimulator battery" on page 98 for more information.

Communication screen

The **Communication** screen appears when the patient programmer is trying to communicate with your neurostimulator. If there is a problem, the displayed screen

Troubleshooting 6

provides instructions (see blue text) on how to resolve the problem and clear the screen.

Table 6.2 describes the **Communication** screen. Unless there is a problem with the communication, the **Communication** screen automatically clears when the process is finished.

Table 6.2 Communication screens

Screen	Description and action
	The patient programmer is communicating with the neurostimulator.
Communication	Continue to hold the patient programmer over your neurostimulator.

Information screens


The information screens provide status information about therapy settings, error conditions, and patient programmer and neurostimulator battery levels.

If the patient programmer audio setting is turned on, three tones alert you to the message.

Table 6.3 describes information screens and instructions on how to proceed (see **blue text**) when these messages appear.

Note: Press any arrow on the **Navigator** key to clear an information screen.

Table 6.3 Information screens

Screen	Description and action
	You tried to increase the amplitude, pulse width, or rate parameters with your therapy off. Turn your therapy on and try communication again.
Turn therapy on	You tried to decrease rate with your therapy off. Turn your therapy on and try decreasing rate again.

Troubleshooting 6

Table 6.3 Information screens (continued)




Screen	Description and action
 Position antenna screen	The patient programmer is trying to communicate with the neurostimulator. Hold the patient programmer or detachable antenna over your neurostimulator.

Table 6.3 Information screens (continued)

Screen	Description and action
 Poor communication	The patient programmer attempted to communicate with the neurostimulator, but communication was unsuccessful. For patients with two neurostimulators, you may have attempted to synchronize with another neurostimulator without first turning the patient programmer off. Reposition the patient programmer over the neurostimulator with the screen facing outward and try communication again. If using the detachable antenna, check that the antenna is connected properly, reposition the antenna and try communication again. Press any arrow on the Navigator key to clear the screen.

Troubleshooting 6

Table 6.3 Information screens (continued)

Screen	Description and action
 Upper limit (amplitude shown)	You tried increasing a parameter (amplitude, pulse width, or rate) above the highest value allowed. Press any arrow on the Navigator key to clear the screen.
 Lower limit (amplitude shown)	You tried decreasing a parameter (amplitude, pulse width, or rate) below the lowest value allowed. Press any arrow on the Navigator key to clear the screen.
 Programmer batteries are low	The patient programmer batteries are low. You can finish programming. Replace the patient programmer batteries before the batteries become depleted.

Troubleshooting 6

Table 6.3 Information screens (continued)





Screen	Description and action
 Power on reset error detected	The neurostimulator has been reset. Your therapy has stopped. Press any arrow on the Navigator key to clear the screen, then press the Therapy on/off  key to turn your therapy on. Call your doctor to report the reset. Note: This information screen has a small information icon in the upper left corner. Compare this POR screen to the warning screens explained on page 110. The warning POR screen has a hazard triangle symbol in the upper left corner.

Table 6.3 Information screens (continued)

Screen	Description and action
 Out of Regulation error	The neurostimulator cannot provide the programmed therapy or perform the requested parameter increase. Reduce the requested parameter increase. (Rechargeable battery only) Recharge your neurostimulator. Call your doctor if these steps do not resolve the issue.
Rechargeable battery only   Rechargeable neurostimulator battery charge level is low	The rechargeable neurostimulator battery charge level is low and therapy will not be available soon. Charge your neurostimulator battery. Refer to the manual packaged with the charging system.

Troubleshooting 6

Table 6.3 Information screens (continued)

Screen	Description and action
Nonrechargeable battery only 	This information screen appears when you synchronize your neurostimulator and patient programmer and the nonrechargeable neurostimulator battery status is low. This means that your therapy will not be available soon. Press any arrow on the Navigator key to clear the screen and call your doctor to report this message.
Nonrechargeable neurostimulator battery status is low	

Possible problems and solutions

Table 6.4 will help you solve problems or identify when to call your doctor. Problems are described in the left column (**bold black text**). The right column lists possible causes of the problem (plain text) and how to correct the problem (**bold blue text**).

Note: If a problem is not solved after several attempts, or if a problem is not described here, contact your doctor.

Table 6.4 Troubleshooting problems

Problems	Causes and actions
Uncomfortable or Intolerable therapy You are experiencing side effects from the therapy.	The selected therapy settings are not suitable for your current activity. 1. Turn your therapy off (see page 44). 2. Call your doctor. If your patient programmer is using Advanced mode, you can try the following before calling your doctor. • Reduce the amplitude or pulse width setting for both sides of your body (see page 66). • If necessary, select a different group (see page 73).

**Table 6.4 Troubleshooting problems
(continued)**

Problems	Causes and actions
Delayed therapy changes You do not notice the therapy right away after turning your therapy on or you notice therapy-like sensations after turning your therapy off.	Your doctor may have programmed SoftStart/Stop™ which gradually starts and stops your therapy. Your doctor can advise you as to how long it takes to feel the effects of the therapy. This time can be different from patient to patient. It also depends upon your medical condition. In general, allow about 8 seconds for your neurostimulator to turn on and off. You may also notice a residual effect after the neurostimulator is turned off.
Intermittent therapy You only notice your therapy some of the time.	Your doctor may have programmed your therapy to turn on and off at regular intervals. However, if you are not receiving adequate symptom relief, contact your doctor.

Troubleshooting 6

**Table 6.4 Troubleshooting problems
(continued)**

Problems	Causes and actions
No therapy You cannot notice your therapy but you think that therapy should be on.	Therapy is off. Use your patient programmer to turn your therapy on (see page 44). The amplitude setting for the active group is set too low to feel. (The effect of the amplitude setting may take up to 10 minutes or more to notice.) If you do not feel the desired effect from a change in amplitude, call your doctor. If your patient programmer is using Advanced mode, you can try the following before calling your doctor. Use your patient programmer to increase the amplitude for one or both sides of your body (see page 69).

**Table 6.4 Troubleshooting problems
(continued)**

Problems	Causes and actions
Patient programmer is unresponsive The display screen is blank when you press a key.	The patient programmer batteries are depleted. Replace the patient programmer batteries (see page 128). The patient programmer batteries are in backwards. Check the battery polarity and reinstall the patient programmer batteries (see page 128).
Dropped patient programmer Your patient programmer falls off a cabinet or table.	The patient programmer is designed to withstand a short drop to a hard surface and still operate normally, even if the case is chipped or nicked. Try the patient programmer; it should work.

**Table 6.4 Troubleshooting problems
(continued)**

Problems	Causes and actions
Fluid on the patient programmer	The patient programmer is not waterproof, and water can damage the device.
Fluid was spilled onto the patient programmer or the patient programmer was dropped into water.	Immediately remove the patient programmer from the water, then dry with a soft towel. Remove the batteries, then allow the battery compartment to air dry at room temperature for 24 hours. Note: If you need to check your rechargeable neurostimulator battery while your patient programmer is drying, you can use the recharger.

Filename Date Time
UC200xxxxxx EN
4.6 x 6 inches (118 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.3: 09-13-2019

7 Maintenance

M929532A059 Rev A 2020-01-01

This chapter describes how to care for your patient programmer and accessories.

Patient programmer batteries

Always keep two new AAA alkaline batteries available for replacement. New batteries provide about two months of use, depending upon how often the patient programmer is used.

You can check your patient programmer batteries at any time. The patient programmer also monitors its own battery charge and will display the following screens when the batteries are low or depleted (Figure 7.1).

△ **Caution:** If the device will not be used for several weeks, remove the batteries from the device. A battery left in the device may corrode, causing damage to the electronic components.

Patient programmer
batteries are low

Patient programmer
batteries are depleted


Figure 7.1 Patient programmer battery screens.

Checking the patient programmer batteries


The images in this procedure are an example
and show how to check the patient
programmer batteries:

1. Synchronize the patient programmer and
neurostimulator.
 - a. Hold the patient programmer (or
antenna) directly over your
neurostimulator with the screen facing
outward.
 - b. Press the **Check** [✓] key. The **Therapy**
screen appears (Figure 7.2).

Maintenance 7

Simple mode

Advanced mode

Figure 7.2 Therapy screen.

2. To display the current patient programmer battery status:

- **Simple mode:** Press the right ► arrow on the **Navigator** key once to display the patient programmer battery status (Figure 7.3).
- **Advanced mode:** Press the up ▲ arrow on the **Navigator** key once to highlight the **Status** row, then press the right ► arrow (Figure 7.3).


Figure 7.3 Status row and Therapy screen.

3. Review the status of the patient programmer battery (Figure 7.4).


Figure 7.4 Patient programmer battery level.

If the batteries are low, see "Replacing the patient programmer batteries" on page 132 for replacement instructions.

Table 7.1 Patient programmer battery replacement screens

	The patient programmer batteries are low. You can finish programming. Press any arrow on the Navigator key to clear the screen; then continue programming. Replace the patient programmer batteries before the batteries become depleted. See page 132.
	The patient programmer batteries are depleted. Programming is not possible. Replace the patient programmer batteries now. See page 132.

Maintenance 7

Replacing the patient programmer batteries

1. Open the battery compartment cover by sliding the cover down and pulling it up.

132 English 37642 2020-01-01

M929532A059 Rev A 2020-01-01


Figure 7.5 Opening the battery cover.

2. Remove the depleted batteries.
3. Insert the new batteries as shown on the battery compartment label.
4. Close the battery compartment cover.
5. Dispose of old batteries according to local requirements.

Maintenance 7

Cleaning and care

Follow these guidelines to ensure that the patient programmer and accessories function properly.

△ **Caution:** If the device will not be used for several weeks, remove the batteries from the device. A battery left in the device may corrode, causing damage to the electronic components.

- Keep the device out of the reach of children and pets. Keep the batteries away from children. If children or pets swallow the batteries, contact a doctor at once.
- Use the device only as explained to you by your doctor or as discussed in this manual.
- Handle the device with care. Do not drop, strike, or step on the device.
- Do not dismantle or tamper with the device.

Maintenance 7

134 English 37642 2020-01-01

M929532A059 Rev A 2020-01-01

- Clean the outside of the device with a damp cloth when necessary. Mild household cleaners will not damage the device or labels.
- The device is not waterproof. Do not allow moisture to get inside the device.
- Keep fresh batteries available.
- Replace low or depleted batteries.

Safety and technical checks

Periodic safety and technical checks or periodic maintenance of the patient programmer are not required. The patient programmer contains no user-serviceable parts. If repair or service is needed, contact your doctor or a Medtronic sales office. Refer to the Medtronic contacts at the end of this manual.

Maintenance 7

Battery and patient programmer disposal

Dispose of depleted batteries and worn out devices according to local requirements. If you no longer need your patient programmer and would like to donate it, contact your doctor.

Declaration of conformity

Medtronic declares that this product is in conformity with the essential requirements of Directive 90/385/EEC on Active Implantable Medical Devices.

For additional information, contact Medtronic. Refer to the list of Medtronic contacts at the end of this manual.

Specifications

Maintenance 7

Table 7.2 Patient programmer specifications

Item	Specification
Power source	2 AAA alkaline batteries (non-rechargeable, LR03)

136 English 37642 2020-01-01

M929532A059 Rev A 2020-01-01

**Table 7.2 Patient programmer specifications
(continued)**

Item	Specification
Operating temperature	9 °C to 43 °C (49 °F to 110 °F)
Temperature limitation ^a	-34 °C to 57 °C (-30 °F to 135 °F)
Size	Approximately 9.4 cm x 5.6 cm x 2.8 cm (3.7 in x 2.2 in x 1.1 in)
Weight, including batteries	Approximately 111 g (3.9 oz.)
Battery life	2 months (average) for alkaline batteries
Service life	Up to 5 years
Mode of operation	Continuous

^a Batteries should be removed from the device for storage or transport.

Maintenance 7

Filename Date Time
UC20000000x EN
4.6 x 6 inches (118 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.3; 09-13-2019

Maintenance 7

138 English 37642 2020-01-01

M929532A059 Rev A 2020-01-01

Filename Date Time
UC200xxxxxx EN
4,6 x 6 inches (118 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.3: 08-13-2019

8 Assistance for the patient programmer

M929532A059 Rev A 2020-01-01

The patient programmer has been designed and tested to provide trouble-free service. If repair or service is needed, contact your doctor or a Medtronic sales office. Refer to the Medtronic contacts at the end of this manual.

[USA] For assistance in the US, call Medtronic patient services at 1-800-510-6735.

If you contact Medtronic about your patient programmer, refer to the serial number. The serial number is located in the battery compartment. This number identifies each patient programmer.

- If your patient programmer stops working, first try the steps in Chapter 6 "Troubleshooting" on page 109. Contact your doctor if indicated by the troubleshooting information. If the problem persists, contact your doctor.
- If you lose your patient programmer, contact your doctor to order a new patient programmer.

Assistance for the patient programmer 8

Filename Date Time
UC200xxxxxx EN
4.6 x 6 inches (118 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.3: 09-13-2019

[USA] Note: To register the patient programmer for service covered by the warranty, complete and mail the warranty registration.

Assistance for the patient programmer 8

37642 2020-01-01 English 141

M929532A059 Rev A 2020-01-01

Index

About rechargeable batteries 46
Adjusting therapy
 adjustment limits 81
 guidelines for 81
 settings 66-82
Adjustment stimulation
 see Adjusting therapy
Alerts
 patient programmer 57, 108
Amplitude
 increasing or decreasing 81, 82
Audio
 tones 66, 111, 115
Backlight (Power/Backlight On/Off) key 34
Batteries (patient programmer)
 battery charge level (low) screen 118
 checking 128
 disposing of 136
 replacing 132
Battery icons 40

Index

142 English 37642 2020-01-01

M929532A059 Rev A 2020-01-01

Battery neurostimulator
 see Neurostimulator battery
Button descriptions 32
Changing a group 73
Charging the Neurostimulator battery
 (rechargeable) 106
Checking batteries
 patient programmer 128
Checking the nonrechargeable battery 47
Communication screen 114
Decreasing therapy settings 69
 limits 81
Disposing of
 patient programmer batteries 136
 patient programmer 136
EOS (end of service) error code 112
Error codes 112
Group, changing 73
Icons
 battery 40
 Parameter/Group 40
 Status 40

Index

37642 2020-01-01 English 143

Identification label 60
Increasing therapy settings
limits 82
Indications 21
Information screen 114
Keys (patient programmer) 32
Low battery level screen
neurostimulator 102
patient programmer 118
Lower limit screens 81
Magnetic resonance imaging (MRI) 84
MRI examinations
activating a new group 86, 91, 92
MRI scan eligibility sheet 85, 88, 93
preparing for an MRI appointment 85
returning therapy to your original group
settings 94
turning off therapy 86, 87, 88
turning therapy back on 89
MRI scan eligibility
types of 84
Navigator key 36

Index

144 English 37642 2020-01-01

M929532A059 Rev A 2020-01-01

Neurostimulator battery 46
Neurostimulator battery (rechargeable)
 automatic alert 57, 108
 charging the battery 106
 checking battery level (charge level) 98
 low battery charge level
 screen 49, 120, 121
 overdischarge 105
 rechargeable battery screens 106
 recharge warning screen 113
 warning 98
 when to check and charge 100
Nonrechargeable neurostimulator battery 47
Options icon 37
Overdischarge 105
Parameter/Group row icons 40
Parameter
 increasing or decreasing limits 81
 limit screens 81
Parameter row 38
Patient programmer alert 57, 108

Index

Patient programmer

batteries 128-133

buttons 32

cleaning and care 134

description 26

disposing of 136

keys 32

low battery level screen 118

purpose 20

replacing batteries 132

screens 110

troubleshooting 110-126

turning on and off 34

POR 96

Possible battery status conditions
(nonrechargeable) 48

Power/Backlight On/Off key 34

Precautions 21

Preference icons 51

Pulse Width

increasing or decreasing 81, 82

Index

146 English 37642 2020-01-01

M929532A059 Rev A 2020-01-01

increasing or decreasing 81, 82

Recharge warning screen 113

Repair 135

Serial number 140

Settings

adjusting therapy 66-82

limits 81

SoftStart/Stop 123

Specifications 136

Status row icons 40

Stimulation

limits when adjusting settings 81

stops suddenly 98

troubleshooting 124

uncomfortable 122

Synchronize

Patient programmer and implanted neurostimulator 41

warning screen 111

37642 2020-01-01 English 147

Therapy settings 37
 decreasing 69
 increasing 69
 turning therapy on or off 33
Therapy
 adjusting settings 66-82
 guidelines for adjusting 81
 On/Off key 33
Tones, audio 66, 111, 115
Turning neurostimulator on or off
 see Turning your therapy on or off
Turning patient programmer on and off 34
Turning your therapy on or off 33
Uncomfortable stimulation
 troubleshooting 122
Upper limit screens 82
Warning screens 111-112
When to check your rechargeable
 neurostimulator battery 100

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