Medtronic

Medtronic

710 Medtronic Parkway
Minneapolis MN 55432-5604
USA

ATTESTATION STATEMENT

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M962281A135 Rev A 2022-06-15

A71100

Restore™ Clinician Programmer Application
Programming Guide

Neurostimulation systems for pain therapy

Neurostimulator models 37022, 37701, 37702, 37703, 37704, 37711, 37712, 37713, 37714, 97702, 97712, 97713, 97714

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
Medtronic

Explanation of symbols

	Conformité Européenne (European Conformity)
	Authorized representative in the European Community
	Manufacturer
	For USA audiences only

Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. AdaptiveStim™, GroupAdjust™, Itrel™, N'Vision™, PrimeAdvanced™, Restore™, RestoreAdvanced™, RestorePrime™, RestoreSensor™, RestoreUltra™, SoftStart/Stop™, SureScan™, and TargetMyStim™ are trademarks of a Medtronic company.

Android™ is a trademark of Google LLC.

The BLUETOOTH® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc., and any use of such marks by Medtronic is under license.

Wi-Fi® and Wi-Fi Alliance® are registered trademarks of Wi-Fi Alliance.

Refer to the indications sheet for indications and related information.

Refer to the device implant manual for device description, package contents, device specifications, and instructions for use.

Refer to the appropriate information for prescribers booklet for contraindications, warnings, cautions, adverse events summary, individualization of treatment, patient selection, use in specific populations, resterilization, and component disposal.

Refer to the MRI Guidelines for Medtronic Neurostimulation Systems for Chronic Pain instructions for use manual for the MRI conditions and MRI-specific warnings and cautions for conducting an MRI scan.

Refer to the System Eligibility, Battery Longevity, Specifications reference manual for neurostimulator selection, battery longevity calculations, and specific neurostimulator specifications.

Refer to the Model 8880T2 Communicator technical manual for warnings, cautions, device description, package contents, device specifications, instructions for use, maintenance information, and the electromagnetic compatibility (EMC) declaration.

Refer to the literature provided by the clinician tablet manufacturer for information regarding wireless use.

 Refer to the clinical summary booklet for information on the clinical study results of the neurostimulation system and individualization of treatment.

Table of Contents

Description of the Programming System.....	7	Setting up the Programming System	12
Model A71100 Restore Clinician Programmer Application ..	7	Preparing the clinician tablet	12
Other components of the programming system	7	Finding and opening the app	12
Programmable Settings	8	Pairing the communicator to the clinician tablet	12
Stimulation parameters.....	8	Positioning the communicator over the neurostimulator....	13
Groups and programs.....	8	Initiating communication with the neurostimulator	13
Optional features	8	Overview of the Clinician Programmer App	15
Data Security and Network Connectivity	10	Using Demo mode.....	15
Data security.....	10	Navigating through workflows.....	15
Network connectivity.....	10	Using the Side Menu	15
Installing application updates	10	On-screen help	16
Installing communicator updates	10	Buttons and indicators on the action bar	16
Returning the clinician tablet	10	When the app is placed in the background	17
General Warnings and Cautions	11	Overview of the Communicator	18
MRI warning for scan eligibility	11	Using the USB connector cable.....	18
MRI warning for trial stimulation	11	Checking battery status of the communicator	18
Sterile field warning for programming components	11	Communicator icons on the tablet status bar	19
EMI caution for telemetry signal disruption.....	11	External Neurostimulator Workflows	20
Cautions related to clinician programming and possible interactions with other devices	11	New ENS: Start Usage > Start Evaluation workflow.....	20
		Configured ENS: Followup workflow	22
		Configured ENS: End Evaluation workflow	23

Implantable Neurostimulator Workflows.....	24	Positioning Leads Relative to Implant Location.....	34
New INS: Start Usage > Implant Device workflow.....	24	Customizing SoftStart/Stop and Cycling	35
Configured INS: Followup workflow.....	26	Changing SoftStart/Stop settings	35
Configured INS: View MRI workflow.....	27	Turning cycling on or off	35
Working with Device and Patient Information ..	28	Enabling Patient Access	36
Checking battery status of the neurostimulator	28	Limits	36
Entering implant location	28	TargetMyStim.....	37
Entering patient information.....	28	Adjustment.....	37
Entering lead information.....	28	Programming Scheduled Therapy	38
Testing lead connectivity.....	29	Checking System Performance.....	39
Entering extension, abandoned lead, or tip location information.....	29	Testing electrode impedance.....	39
Clearing information in external neurostimulator for re-use	29	Troubleshooting electrode impedance results	40
Working with Groups and Programs	30	Testing group impedance	41
Programs screen	30	Troubleshooting group impedance results.....	41
Program subscreen	30	Preparing System for an MRI Scan.....	43
Turning stimulation on or off	30	View MRI workflow	43
Adjusting Stimulation Parameters.....	31	AdaptiveStim Technology.....	44
High-output interlocks.....	31	Caution regarding cervical location	44
Shared electrodes	31	Caution regarding healing period	44
Assigning electrode configuration.....	31	Configuring AdaptiveStim Technology with the wizard	44
Assigning pulse width, rate, and amplitude	32		

Reconfiguring an AdaptiveStim-enabled group	45	Working with Reports	55
What you can do on the AdaptiveStim screen	45	Using the Reports screen	55
Reviewing position amplitudes for all AdaptiveStim-enabled groups	48	Patient Data Service app.....	56
Testing AdaptiveStim settings	48	Finding and opening the app	56
Adjusting position amplitudes for individual positions	48	Accessing or deleting a patient session record	56
Checking patient positions	49	Deleting multiple patient session records at once	56
Reorienting patient positions	49	Setting the PDF auto-delete feature	57
Disabling AdaptiveStim Technology for a group	49	Viewing System Information	58
Enabling or disabling AdaptiveStim Technology for the patient	49	Viewing information about system components	58
Turning the Position Trend diary on or off	50	Viewing information about trademarks and licenses	58
Adjusting position angles	50	Changing device date/time or implant date information	58
Adjusting transition times	50	Compiling a Medtronic Data Report	58
Adjusting mobility rate	51	Viewing Alerts	59
Turning off mobility detection	51	Working with alerts	59
Troubleshooting AdaptiveStim Technology	51	Troubleshooting	60
Viewing Diaries	54	Data validation error message	60
Stimulation Usage	54	Power on reset (POR) message	60
Group Usage	54	Out of regulation (OOR) message	60
Position Trend	54	Overcharge message	60
Adjustments	54	Overdischarge message	60
Resting Trend	54	Rapid discharge message	60
Recharge	54	Elective Replacement Indicator (ERI) message	60

End Of Service (EOS) message.....	60
Using impedance measurements for troubleshooting	61
Communicator and clinician tablet troubleshooting	61
Physician Recharge mode.....	64
Patient programmer messages	65

Description of the Programming System

Model A71100 Restore Clinician Programmer Application

The Medtronic Model A71100 Restore Clinician Programmer Application (app) is intended for use by clinicians in the programming of the following Medtronic neurostimulators for pain therapy.

External Neurostimulator (ENS):

- Model 37022 External Neurostimulator

Implantable Neurostimulators (INS):

- Model 37703 Itrel 4 Implantable Neurostimulator
- Model 37704 Itrel 4 Implantable Neurostimulator
- Model 37711 Restore Implantable Neurostimulator
- Model 37701 RestorePrime Implantable Neurostimulator
- Model 37702 PrimeAdvanced Implantable Neurostimulator
- Model 97702 PrimeAdvanced SureScan MRI Implantable Neurostimulator
- Model 37712 RestoreUltra Implantable Neurostimulator
- Model 97712 RestoreUltra SureScan MRI Implantable Neurostimulator
- Model 37713 RestoreAdvanced Implantable Neurostimulator
- Model 97713 RestoreAdvanced SureScan MRI Implantable Neurostimulator
- Model 37714 RestoreSensor Implantable Neurostimulator
- Model 97714 RestoreSensor SureScan MRI Implantable Neurostimulator

Other components of the programming system

The clinician programmer app is intended for use with the following components.

Model CT900 Clinician Tablet

The clinician programmer app is loaded on the Model CT900 Clinician Tablet with Android-based operating system.

Model 8880T2 Communicator

The Model 8880T2 Communicator is a handheld device used in conjunction with the clinician tablet and clinician programmer app to communicate with Medtronic neurostimulators.

Model A901 Communication Manager Application

The Model A901 Communication Manager Application is an application on the clinician tablet that manages the telemetry communications for the clinician tablet.

Model A902 Patient Data Service Application

The Model A902 Patient Data Service Application is an application on the clinician tablet that manages the storage of patient session data on the clinician tablet. The Patient Data Service app can be used to access or delete patient session data for any patients whose Medtronic devices have been programmed using any clinician programmer app on the clinician tablet.

Programmable Settings

The clinician programmer app is used to enter, review, and modify programmable settings in the neurostimulator. These programmable settings define the stimulation therapy delivered to the patient.

Stimulation parameters

Stimulation is delivered in electrical pulses. Stimulation parameters define the attributes of these pulses and can be adjusted to manage patient therapy.

- **Electrode configuration:** Programmed negative (–) polarity, positive (+) polarity, or off.
- **Pulse width:** Duration of each pulse in microseconds (µs).
- **Rate:** Frequency of pulses in hertz (Hz).
- **Amplitude:** Strength of a pulse in volts (V).

Groups and programs

A program is a specific combination of pulse width, rate, and amplitude settings acting on a specific electrode configuration. The program defines the stimulation pulses that will be delivered for therapy. For certain neurostimulator models, multiple programs can be combined into groups. Groups are designated as A, B, C, and so on in the programmer app, while programs are designated as 1, 2, 3, 4. Refer to Table 1 on page 9 for a summary of the group and program options available for each neurostimulator model.

When using more than one program, the pulses are delivered sequentially—first a pulse from one program, then a pulse from the next program, and so on. Pulses from different programs are never delivered simultaneously. Each group and its associated programs can be used to provide therapy for specific areas of coverage or specific patient activities.

Pulse width, amplitude, and electrode configuration are assigned at the program level: that is, each program within a group can have different values. Rate is assigned at the group level: that is, each program within a group will have the same rate.

Optional features

Optional features are available for tailoring the stimulation to the patient's needs. Refer to Table 1 on page 9 for a summary of the features supported in each neurostimulator model.

- **Cycling:** Cycles stimulation on and off at clinician-determined intervals.
- **SoftStart/Stop:** Slowly increases the amplitude when stimulation is turned on and slowly decreases the amplitude when stimulation is turned off.
- **Patient Limits:** Identifies the stimulation parameters that the patient can adjust and sets the upper limits for those parameters.
- **Adjustment (also known as GroupAdjust):** Allows the patient to increase or decrease all program amplitudes in the active group at the same time.
- **TargetMyStim:** Allows the patient to move stimulation settings from electrode to electrode up and down the lead.
- **Scheduled Therapy:** Schedules specific groups to be active at programmed times during a 24-hour period.
- **AdaptiveStim Technology:** Adjusts stimulation automatically when the patient changes position.

Table 1. Neurostimulator features.

Model	Leads	Electrodes	Groups	Programs/Group	Total Programs	Cycling	SoftStart/Stop	Patient Limits	Adjustment	TargetMyStim	Scheduled Therapy	AdaptiveStim
37022 ENS	4	16	8	4	16	●	●	●	●	●	●	–
37703 / 37704 Itrel 4	1	4	1	1	1	●	●	●	–	–	–	–
37711 Restore	2	16	26	4	32	●	●	●	–	–	●	–
37701 RestorePrime	2	16	26	4	32	●	●	●	–	–	●	–
37702 PrimeAdvanced	4	16	26	4	32	●	●	●	●	●	●	–
97702 PrimeAdvancedMRI	4	16	26	4	32	●	●	●	●	●	●	–
37712 RestoreUltra	4	16	8	4	16	●	●	●	●	●	●	–
97712 RestoreUltraMRI	4	16	8	4	16	●	●	●	●	●	●	–
37713 RestoreAdvanced	4	16	26	4	32	●	●	●	●	●	●	–
97713 RestoreAdvancedMRI	4	16	26	4	32	●	●	●	●	●	●	–
37714 RestoreSensor	4	16	8	4	16	●	●	●	●	●	–	●
97714 RestoreSensorMRI	4	16	8	4	16	●	●	●	●	●	–	●

Legend: ● Feature is available. – Feature is not available.

Data Security and Network Connectivity

Data security

The clinician programming system uses and stores data about the patient's health and implanted medical device. This data is protected by application-level encryption and encryption provided by the clinician tablet. The clinician programming system does not provide data protection for data exported to another destination. Exported data should be handled in accordance with your facility's security policy for data handling and storage.

- Medtronic recommends that you always save exported data to the default reports location on the clinician tablet.

Network connectivity

Network connectivity is required for initial app registration and for installation of Medtronic app updates and communicator firmware updates. Network connectivity is not required for neurostimulator programming. To protect your clinician programming system, Medtronic recommends you implement the following security measures:

- Secure your clinician tablet by disabling network connectivity during any programming session.
- Use a managed, trusted Wi-Fi connection when network connectivity is needed.
- Connect the clinician tablet to the network periodically to check for update notifications.

 **Caution:** Connecting the clinician tablet to a network that includes other equipment could result in unforeseen risks to patients, operators, or third parties. Changes to your network such as adding,

disconnecting, and upgrading equipment; upgrading or installing software; or changing network configurations could also introduce additional risks. Analyze, evaluate, and control any identified risks.

If you suspect a cybersecurity event has occurred, stop using the app (if possible) and contact your IT Security team or Medtronic Technical Services to document and respond to the suspected incident.

If your clinician tablet is lost or stolen, contact Medtronic using the contact information listed on the back cover of this manual.

Installing application updates

Medtronic periodically updates the therapy app and will not install updates without notifying you. Network connectivity is required to update the app. When notified that an app update is available, follow the instructions provided by Medtronic to install the update.

Installing communicator updates

When notified to install a communicator firmware update, connect the communicator to the clinician tablet using the USB connector cable. Follow the instructions provided by Medtronic to install the update.

Returning the clinician tablet

If you need to return the clinician tablet for disposal or replacement, contact Medtronic using the contact information listed on the back cover of this manual. Instructions will be provided for preparing the clinician tablet for return.

General Warnings and Cautions

If a serious incident related to the app occurs, immediately report the incident to Medtronic and the applicable competent authority.

MRI warning for scan eligibility

 **Warning:** For neurostimulators with SureScan MRI Technology (97702, 97712, 97713, 97714), always program the following component information:

- Neurostimulator and lead model numbers
- Neurostimulator implant location and lead tip location
- Presence of extensions or abandoned leads

If this information is not up-to-date or is entered incorrectly, MRI scan eligibility data will be inaccurate, and the patient is at risk for one of the following:

- The patient is allowed to have an MRI scan inappropriate for the implanted components, which could cause tissue heating, resulting in tissue damage or serious patient injury.
- The patient is unnecessarily restricted from having an MRI scan.

For all neurostimulators, refer to the *MRI Guidelines for Medtronic Neurostimulation Systems for Chronic Pain instructions for use* manual for MRI scan eligibility information, MRI scan conditions, and MRI-specific warnings and cautions for conducting an MRI scan.

MRI warning for trial stimulation

 **Warning:** Physicians should not prescribe MRI for patients undergoing trial stimulation or who have any neurostimulation system components that are not fully implanted. Explant all trial stimulation components if an MRI scan is required. MRI has not been tested

on trial stimulation components and may cause heating of the lead electrodes, resulting in tissue damage or serious patient injury.

Sterile field warning for programming components

 **Warning:** To use the programming components (external and nonsterile) in a sterile field, place a sterile barrier between the patient and the programming components to prevent infection. Do not sterilize any of the programming components. Sterilization may damage the programming components.

EMI caution for telemetry signal disruption

 **Caution:** Electromagnetic interference (EMI) can disrupt programming and telemetry communication. If EMI interference is suspected, such as from radio frequency identification (RFID) equipment, move away from the likely source of EMI and try again.

Cautions related to clinician programming and possible interactions with other devices

Refer to the *Information for Prescribers* booklet for cautions related to clinician programming and possible interactions with other devices (such as cardiac devices).

Setting up the Programming System


Figure 1. Components of the programming system.

Preparing the clinician tablet

Refer to the *Model CT900 Clinician Tablet quick start guide* for instructions on how to turn on the tablet and complete initial setup.

Check the battery level of the clinician tablet.

- The tablet's battery level is shown on the tablet status bar (uppermost row on the tablet screen).
- Recharge the tablet if the battery level is low.

Finding and opening the app

1. Navigate to the Apps on the clinician tablet.
2. Find the **Restore** app icon.
3. Tap the **Restore** app icon to open the app.

Pairing the communicator to the clinician tablet

Refer to the *Model 8880T2 Communicator technical manual* for additional instructions, if needed. Only one communicator can be paired with a clinician tablet at a time.

1. Connect the USB connector cable to the communicator and the clinician tablet.
2. Turn on the communicator (slide the power button down, then release).
3. A message appears on the clinician tablet, which asks permission to open the Communication Manager when the USB device is connected:
 - a. Select to use the USB device by default. This selection prevents this message from appearing again.
 - b. Select OK.
4. On the initial app screen, tap **CONNECT**. When pairing is complete, the LED light (■ — ■) on the communicator turns a solid green.

Once paired, you can disconnect the USB connector cable. The communicator and clinician tablet communicate using BLUETOOTH® wireless technology. The clinician tablet must remain within 3 meters (10 feet) of the communicator for wireless communication.

Positioning the communicator over the neurostimulator

The communicator uses near-field magnetic induction communication (also referred to as "proximal telemetry") to communicate with the external and implantable neurostimulators. This telemetry requires close proximity for successful programming.

1. Hold the communicator directly over the neurostimulator, and position the target symbol (⊕) on the back of the communicator so that it is centered over and facing the neurostimulator (Figure 2).
 - The target symbol on the communicator indicates the location of the internal antenna.


Figure 2. Position target symbol on communicator over neurostimulator (INS shown, patient supine).

Note: For the external neurostimulator, you can either hold the communicator over the neurostimulator as shown in Figure 2 or you

can attach the neurostimulator directly to the communicator as shown in Figure 3.

- Locate the slot on the front of the communicator. Place the external neurostimulator into the slot, and rotate it so the output jack faces toward the outer edge of the communicator. Then attach the trialing cable.


Figure 3. External neurostimulator attached to communicator.

Note: Avoid placing the communicator directly on a metal surface. Metal surfaces can interfere with communication between the communicator and the neurostimulator.

Initiating communication with the neurostimulator

1. On the Search for Device screen, tap FIND DEVICE.
Or
Press the Communicate button (↔) on the communicator.
2. Select the serial number of the neurostimulator you intend to configure, and tap CONNECT.

3. After connecting to the neurostimulator, the **CURRENT DEVICE STATUS** screen will display.

From this screen, you can select workflow options. Refer to the following sections for workflow details:

- "External Neurostimulator Workflows" on page 20.
- "Implantable Neurostimulator Workflows" on page 24.

Notes:

- If you have problems establishing communication, refer to "Communicator and clinician tablet troubleshooting" on page 61.

- Do not leave the clinician tablet unattended during an active programming session.
- When a configured neurostimulator is interrogated, a lead connectivity test is performed to check for impedance issues. The patient may feel a stimulation change during the test.
- If the last programming session was completed with the Model 8840 N'Vision Clinician Programmer, programming changes may occur upon interrogation with the app. Refer to Table 2 for a list of the potential programming changes.

Table 2. Potential programming changes upon interrogation (if last programmed by N'Vision Clinician Programmer).

Feature	Programming Change	Notes
Amplitude Resolution	Amplitude resolution between 0 V and 10 V will be set to 0.05 V for all programs.	The patient will have to press the Increase or Decrease key twice for each change of 0.1 V.
Patient Limits	All programmed patient limits will be set to Custom limit type. Lower patient limits will be set to the following minimum values: Amplitude 0.0 V, Pulse Width 60 μ s, Rate 2 Hz.	The app only supports Custom limit type and upper limit customization for patient limits. Existing values for upper patient limits will not change upon interrogation.
SoftStart/Stop	All groups will be set to the highest SoftStart/Stop value previously configured in the neurostimulator.	The app supports a single SoftStart/Stop value for all groups. If some groups are On and some groups are Off, the setting will change to On for all groups.
Cycling On/Off Times	Cycling On/Off times will be set to the highest programmed values for all groups with Cycling enabled.	The app supports a single set of Cycling On/Off times for all groups. Cycling can be enabled/disabled at the group level.
AdaptiveStim Technology	For groups with AdaptiveStim enabled, all positions (except Mobile) will be enabled. If Mobile was previously enabled, there will be no change to that position. The position amplitude for newly enabled positions will be set to 0.0 V.	The patient may perceive therapy changes due to the newly enabled positions. Make sure to review all AdaptiveStim settings before ending the programming session.

Note: Alert messages in the app will provide details of the specific programming changes to the neurostimulator.

Overview of the Clinician Programmer App

Using Demo mode

Use Demo mode to explore the clinician programmer app without interrogating or updating a neurostimulator.

When the app is in Demo mode, it is not communicating with a real neurostimulator and does not require the use of a communicator. Any data displayed in Demo mode are not actual.

Demo mode can be used for training and demonstration purposes and to familiarize yourself with the app interface before starting an actual programming session.

1. On the initial app screen, tap **DEMO**.
2. Select a device, and tap **CONNECT**.
3. Navigate through the workflow. One workflow at a time can be explored in a Demo session.
4. To exit Demo mode, go to the **Summary** screen and tap **EXIT WORKFLOW**.

Exiting the workflow will return you to the initial app screen where you can select **DEMO** again and select another device to explore.

Navigating through workflows

The screens typically used during a programming session are grouped into workflows. The workflows vary depending on the device type (external neurostimulator vs. implantable neurostimulator) and device status (new vs. configured). The workflows may also vary depending on neurostimulator features, such as SureScan MRI Technology.

After you start a workflow, a workflow navigator appears at the top of the screen. For example, the Configured ENS: Followup workflow contains these screens:

Device ► Lead Select ► Tip Location ► Programs ► Electrode Impedance ► Diaries ► Reports ► Summary

- To move from screen to screen, swipe left or right on the screen.
- Alternately, you can tap on the screen name in the workflow navigator to go directly to that screen.
 - If the workflow navigator extends beyond the width of the screen, swipe the workflow navigator left or right to show the screens that are clipped off.
- The current screen is underlined and highlighted.
- If input is required in a workflow screen, you will not be able to advance to the next screen until you enter the required input.
 - The app will prompt you for the required input.
 - Red asterisks indicate required fields.

Using the Side Menu

Screens not in the current workflow can be accessed from the Side Menu. For example, the **Patient Info** screen is in the new device workflow, but not in the configured device workflow. To change patient data in a configured device, you would access the **Patient Info** screen from the Side Menu.

- To access the Side Menu, tap the **Side Menu** button (≡) in the top left corner on the action bar.
- Alternately, you can swipe from the left edge of the tablet to access the Side Menu.

On-screen help

Help information is available in the clinician programmer app for workflow screens. View the help information to learn how the app interface works and how to complete tasks on the screen.

- To display help information, tap the Help button (?) in the top right corner on the action bar.
- Tap the Help button (or tap anywhere on the screen) to close the help screen.

Buttons and indicators on the action bar

An action bar appears at the top of all screens in all workflows. Table 3 describes the buttons and indicators that may appear on the action bar.

Table 3. Buttons and indicators.

Button / Indicator	Description
	Side Menu button: Tap to access screens outside of the workflow. (Alternately, you can swipe from the left edge of the tablet to access the Side Menu.)
	Back button: Tap to close a screen you accessed from the Side Menu and return to the previous workflow screen.
DEMO:	Demo mode indicator: Appears in the action bar when in Demo mode.

Button / Indicator	Description
	Stimulation toggle: Tap to turn stimulation on or off. Green indicates that stimulation is on. Appears in the action bar when at least one program exists.
	AdaptiveStim indicator: Applicable to neurostimulators with AdaptiveStim Technology (37714, 97714). Appears in the action bar when AdaptiveStim Technology is enabled.
	MRI Status indicator: Applicable to neurostimulators with SureScan MRI Technology (97702, 97712, 97713, 97714). Appears in the action bar when stimulation is off and the device is ready for an MRI scan. The MRI Status indicator does not provide MRI scan eligibility information.
	Help button: Tap to display help information for the current screen. Tap the Help button (or tap anywhere on the screen) to close the help screen.
	Go to Summary screen button: Tap to go directly to the Summary screen, and then exit the workflow to properly end the programming session.

When the app is placed in the background

When you navigate away from the clinician programmer app while in a programming session, the app is placed in the background and is represented as a floating widget on the clinician tablet.

The floating widget is an overlay that floats on top of the screen to remind you that you are still in a session.

- Press and drag the floating widget to move it on the tablet screen.
- Tap on the floating widget to access the app and resume the session.

Overview of the Communicator

The Model 8880T2 Communicator is a non-sterile component used in conjunction with the clinician tablet and clinician programmer app to communicate with Medtronic neurostimulators.

 **Warning:** To use the communicator in a sterile field, place a sterile barrier between the patient and the communicator to prevent infection. Do not sterilize any of the programming components. Sterilization may damage the programming components.

The communicator is handheld and battery-operated. Communication between the communicator and the clinician tablet can occur wirelessly using BLUETOOTH® technology or wired using the USB connector cable.

For information on the Model 8880T2 Communicator, including description, specifications, components, instructions for use, maintenance, and troubleshooting, refer to the *Model 8880T2 Communicator technical manual*.

Using the USB connector cable

The following situations require the use of the USB connector cable:

- First-time pairing of any communicator with a clinician tablet.
- Firmware updates to the communicator.
- Environments where multiple devices are using BLUETOOTH® wireless technology and thereby creating interference.
- Environments where BLUETOOTH® wireless technology is prohibited.

Checking battery status of the communicator

Check the battery status of the communicator before a programming session and regularly during a programming session. When you turn on the communicator, a battery LED indicator displays on the top right side of the communicator.

- **Solid Green:** Communicator has been turned on. Battery level is acceptable and all communicator functions are enabled.
- **Solid Amber:** Communicator has been turned on. Batteries should be replaced soon.
- **Flashing Red:** Communicator has been turned on. Batteries should be replaced immediately. **Communicate** button (↔) is disabled.

When the communicator is paired to the clinician tablet, swipe down from the top of the tablet to view the following during a programming session:

- The battery level of the communicator.
- The status of the connection between the communicator and the clinician tablet.

Communicator icons on the tablet status bar

Table 4 describes the icons that appear in the tablet status bar (uppermost row on the tablet screen) indicating the status of the connection between the communicator and the clinician tablet.

Table 4. Communicator icons in the tablet status bar.

Communicator Icon	Description
	The communicator is in the process of connecting to the clinician tablet.
	The communicator is communicating with the clinician tablet via the USB connector cable.
	The communicator requires the USB connector cable to communicate with the clinician tablet.
	The communicator is communicating with the clinician tablet using BLUETOOTH® wireless technology.

External Neurostimulator Workflows

When an external neurostimulator (ENS) is interrogated, the clinician programmer app determines if the ENS is new or configured. An ENS is identified as 'new' if it is in shelf state or if it has been cleared for re-use. An ENS is identified as 'configured' if it has patient data stored in it. Based on the status of the ENS, the clinician programmer app automatically presents the applicable workflow options.

New ENS:

- **Start Usage > Start Evaluation workflow – go to page 20.**

Use at the start of a trial period, when the leads are implanted and the ENS is connected.

Use for testing lead placement.

Configured ENS:

- **Followup workflow – go to page 22.**

Use at followup sessions with the patient during the trial period.

- **End Evaluation workflow – go to page 23.**

Use at the last followup session with the patient at the end of the trial period.

Use to clear the ENS for re-use.

Notes:

- To display on-screen help information, tap the **Help** button (🔗).
- The screens typically used during a programming session are shown in the workflow navigator. Additional screens not in the current workflow can be accessed from the Side Menu (☰).
- Unsaved information may be lost if the application is improperly terminated during a programming session. To end a session, go to the **Summary** screen and tap the **EXIT WORKFLOW**

button. Ending a session generates a complete session record that is stored in a database on the clinician tablet.

New ENS: Start Usage > Start Evaluation workflow

1. On the **CURRENT DEVICE STATUS** screen, confirm that the device information corresponds with the intended ENS.
2. Tap **START USAGE**.
This action sets the In Use Date to the current date.
3. Under **START EVALUATION**, tap **START**.
4. You are now in the **New ENS: Start Evaluation** workflow. Use Table 5 to guide you through the screens and tasks in the workflow.

Table 5. New ENS: Start Evaluation workflow.

Screen	Tasks	More information
Device	Check the battery status of the neurostimulator.	Page 28
Patient Info	Enter patient information.	Page 28
Lead Select	Select the implanted leads. Select a lead from the scrollable drop-down list, and then drag and drop it at the appropriate electrode numbers. Tap UPDATE to save the data.	Page 28
	Ensure leads are properly connected to the neurostimulator, and then check connectivity by tapping the CHECK CONNECTIVITY button.	Page 29
Tip Location	Indicate whether extensions are implanted.	Page 29
Programs (and subscreens)	Create programs in groups. Tap a group to make it active. A highlighted box around a group indicates the group is active. Tap the plus sign to add a program within a group.	Page 30
	Program subscreen: Assign the electrode configuration for the lead(s). Assign pulse width, rate, and amplitude.	Page 30
	Lead Manipulation subscreen (only applicable if more than one lead): Use to represent the relative positioning of the implanted leads as viewed from the patient's back.	Page 34
	Energy subscreen: Set parameters for SoftStart/Stop and cycling.	Page 35
	Patient Access subscreen: Assign the stimulation parameters that the patient can adjust.	Page 36
Electrode Impedance	Measure electrode impedance. These measurements verify the integrity of lead, extension, and connector pathways. They may provide information about lead problems (such as lead breakage, short circuit, open circuit).	Page 39
Group Impedance	Measure group (therapy) impedance. These measurements provide documentation that the pathways are intact and the current provided is sufficient for the selected therapy.	Page 41
Summary	Confirm that the therapy settings shown are what is intended for the patient. End the session properly by tapping the EXIT WORKFLOW button.	--

Configured ENS: Followup workflow

1. On the **CURRENT DEVICE STATUS** screen, confirm that the patient and device information corresponds with the intended patient and ENS.

2. Under **FOLLOWUP**, tap **START**.

3. You are now in the **Configured ENS: Followup** workflow. Use Table 6 to guide you through the screens and tasks in the workflow.

Table 6. Configured ENS: Followup workflow.

Screen	Tasks	More information
Device	Check the battery status of the neurostimulator.	Page 28
Lead Select	Review the information, and make any updates. Updates may delete programming.	Page 28
Tip Location	Review the information, and make any updates.	Page 29
Programs (and subscreens)	Review the groups, programs, and settings, and make any updates. Confirm that the active group is the intended group.	Page 30
	Program subscreen: Assign the electrode configuration for the lead(s). Assign pulse width, rate, and amplitude.	Page 30
	Lead Manipulation subscreen (only applicable if more than one lead): Use to represent the relative positioning of the implanted leads as viewed from the patient's back.	Page 34
	Energy subscreen: Set parameters for SoftStart/Stop and cycling.	Page 35
	Patient Access subscreen: Assign the stimulation parameters that the patient can adjust.	Page 36
Electrode Impedance	Measure electrode impedance. These measurements verify the integrity of lead, extension, and connector pathways. They may provide information about lead problems (such as lead breakage, short circuit, open circuit).	Page 39
Diaries	View the Stimulation Usage and Group Usage diaries to assess the patient's use of the therapy.	Page 54
Reports	Select a report to view, download, or delete.	Page 55
Summary	Confirm that the therapy settings shown are what is intended for the patient. End the session properly by tapping the EXIT WORKFLOW button.	--

Configured ENS: End Evaluation workflow

1. On the **CURRENT DEVICE STATUS** screen, confirm that the patient and device information corresponds with the intended patient and ENS.
2. Tap **END EVALUATION > START**.
3. You are now in the **Configured ENS: End Evaluation** workflow. Use Table 7 to guide you through the screens and tasks in the workflow.

Table 7. Configured ENS: End Evaluation workflow.

Screen	Tasks	More information
Device	Check the battery status of the neurostimulator.	Page 28
Lead Select	Review the information, and make any updates. Updates may delete programming.	Page 28
Tip Location	Review the information, and make any updates.	Page 29
Programs	Review the groups, programs, and settings, and make any updates.	Page 30
Diaries	View the Stimulation Usage and Group Usage diaries to assess the patient's use of the therapy.	Page 54
Reports	Select a report to view, download, or delete.	Page 55
Clear Device	Clear the patient and programming information from the device. This action prepares the external neurostimulator for re-use.	Page 29
Summary	End the session properly by tapping the EXIT WORKFLOW button.	--

Implantable Neurostimulator Workflows

When an implantable neurostimulator (INS) is interrogated, the clinician programmer app determines if the INS is new or if it has been previously configured. Based on the status of the INS, the clinician programmer app automatically presents the applicable workflow options.

New INS:

- **Start Usage > Implant Device workflow – go to page 24.**
Use when a new INS is being configured (programmed) for the first time.

Configured INS:

- **Followup workflow – go to page 26.**
Use at followup sessions with the patient.
- **View MRI workflow – go to page 27.**
The View MRI workflow is only applicable to neurostimulators with SureScan MRI Technology (97702, 97712, 97713, 97714). Use to determine MRI scan eligibility and prepare a patient's neurostimulation system for an MRI scan.

Notes:

- To display on-screen help information, tap the **Help** button (🔗).
- The screens typically used during a programming session are shown in the workflow navigator. Additional screens not in the current workflow can be accessed from the Side Menu (☰).
- Unsaved information may be lost if the application is improperly terminated during a programming session. To end a session, go to the **Summary** screen and tap the **EXIT WORKFLOW** button. Ending a session generates a complete session record that is stored in a database on the clinician tablet.

New INS: Start Usage > Implant Device workflow

You can initiate communication with the INS in the package before it is moved into the sterile field. Tasks such as checking the battery status and entering device and patient information can be performed without removing the INS from the package.

1. On the **CURRENT DEVICE STATUS** screen, confirm that the device information corresponds with the intended INS.

2. Tap **START USAGE**.

This action takes the neurostimulator out of shelf state and sets the In Use Date/Implant Date to the current date. For rechargeable neurostimulators, this action starts a clock inside the neurostimulator that determines the date when the Elective Replacement Indicator (ERI) and End Of Service (EOS) messages appear.

The In Use Date and Implant Date are the same date. The In Use Date is shown on the **CURRENT DEVICE STATUS** screen, and the Implant Date is shown on the **About System** screen. You can change the Implant Date on the **About System** screen, but it will not impact the clock inside the neurostimulator that determines ERI and EOS.

3. Under **IMPLANT DEVICE**, tap **START**.
4. You are now in the **New INS: Implant Device** workflow. Use Table 8 to guide you through the screens and tasks in the workflow.

Table 8. New INS: Implant Device workflow.

Screen	Tasks	More Information
Device	Check the battery status of the neurostimulator.	Page 28
	Select the implant location.	Page 28
Patient Info	Enter patient information.	Page 28
Lead Select	Select the implanted leads. Select a lead from the scrollable drop-down list, and then drag and drop it at the appropriate electrode numbers. Tap UPDATE to save the data.	Page 28
	Ensure leads are properly connected to the neurostimulator, and then check connectivity by tapping the CHECK CONNECTIVITY button.	Page 29
Tip Location	Indicate whether extensions are implanted. Indicate whether abandoned leads remain implanted (for MRI models only). Select the tip location for each lead (for MRI models only). (This screen is not supported for Itrel 4.)	Page 29
Programs (and subscreens)	Create programs in groups. Tap a group to make it active. A highlighted box around a group indicates the group is active. Tap the plus sign to add a program within a group.	Page 30
	Program subscreen: Assign the electrode configuration for the lead(s). Assign pulse width, rate, and amplitude.	Page 30
	Lead Manipulation subscreen (only applicable if more than one lead): Use to represent the relative positioning of the implanted leads as viewed from the patient's back.	Page 34
	Energy subscreen: Set parameters for SoftStart/Stop and cycling.	Page 35
	Patient Access subscreen: Assign the stimulation parameters that the patient can adjust.	Page 36
Electrode Impedance	Measure electrode impedance. These measurements verify the integrity of lead, extension, and connector pathways. They may provide information about lead problems (such as lead breakage, short circuit, open circuit).	Page 39
Group Impedance	Measure group (therapy) impedance. These measurements provide documentation that the pathways are intact and the current provided is sufficient for the selected therapy.	Page 41
Summary	Confirm that the therapy settings shown are what is intended for the patient. End the session properly by tapping the EXIT WORKFLOW button.	--

Configured INS: Followup workflow

1. On the **CURRENT DEVICE STATUS** screen, confirm that the patient and device information corresponds with the intended patient and INS.

2. Under **FOLLOWUP**, tap **START**.

3. You are now in the **Configured INS: Followup** workflow. Use Table 9 to guide you through the screens and tasks in the workflow.

Table 9. Configured INS: Followup workflow.

Screen	Tasks	More information
Device	Check the battery status of the neurostimulator.	Page 28
Lead Select	Review the information, and make any updates. Updates may delete programming.	Page 28
Tip Location	Review the information, and make any updates. (This screen is not supported for Itrel 4.)	Page 29
Programs (and subscreens)	Review the groups, programs, and settings, and make any updates. Confirm that the active group is the intended group.	Page 30
	Program subscreen: Assign the electrode configuration for the lead(s). Assign pulse width, rate, and amplitude.	Page 30
	Lead Manipulation subscreen (only applicable if more than one lead): Use to represent the relative positioning of the implanted leads as viewed from the patient's back.	Page 34
	Energy subscreen: Set parameters for SoftStart/Stop and cycling.	Page 35
	Patient Access subscreen: Assign the stimulation parameters that the patient can adjust.	Page 36
Electrode Impedance	Measure electrode impedance. These measurements verify the integrity of lead, extension, and connector pathways. They may provide information about lead problems (such as lead breakage, short circuit, open circuit).	Page 39
AdaptiveStim	Configure the settings for AdaptiveStim Technology (RestoreSensor models 37714, 97714 only).	Page 44
Diaries	Select a diary to view from the drop-down list. Diaries provide information on therapy usage. (This screen is not supported for Itrel 4.)	Page 54
Reports	Select a report to view, download, or delete.	Page 55
Summary	Confirm that the therapy settings shown are what is intended for the patient. End the session properly by tapping the EXIT WORKFLOW button.	--

Configured INS: View MRI workflow

The View MRI workflow is only applicable to neurostimulators with SureScan MRI Technology (97702, 97712, 97713, 97714).

1. On the **CURRENT DEVICE STATUS** screen, confirm that the patient and device information corresponds with the intended patient and INS.
2. Tap **VIEW MRI > START**.
3. You are now in the **Configured INS: View MRI** workflow. Use Table 10 to guide you through the screens and tasks in the workflow.

Table 10. Configured INS: View MRI workflow.

Screen	Tasks	More information
Device	Confirm the neurostimulator information and implant location.	Page 43
Patient Info	Confirm the patient information.	Page 43
Lead Select	Confirm the implanted lead information.	Page 43
Tip Location	Confirm the extension, abandoned lead, and tip location information.	Page 43
MRI	Review the MRI scan eligibility. Determine whether the patient is eligible for the intended MRI scan. If patient is eligible, turn stimulation off to prepare the patient's neurostimulation system for the scan.	Page 43
Reports	Prepare the MRI Report. The MRI Report contains information intended for MRI clinicians about the MRI scan eligibility of the implanted system.	Page 43
Summary	Confirm that stimulation is off and MRI status shows the device is ready for an MRI scan. End the session properly by tapping the EXIT WORKFLOW button.	--

Working with Device and Patient Information

Checking battery status of the neurostimulator

The battery status for the neurostimulator (ENS or INS) is shown at the start and at the end of a programming session on the following screens:

- CURRENT DEVICE STATUS screen
- Device screen
- Summary screen

For all neurostimulators, the battery status is displayed in volts. For rechargeable neurostimulators, a battery image also displays representing the battery level percentage: 0%, 25%, 50%, 75%, or 100%.

At the start of a programming session, if the battery voltage is below 25% for a rechargeable neurostimulator or an ENS, a message notifies you that the device is running on low battery.

If the battery voltage for a rechargeable neurostimulator or an ENS is too low to sustain therapy or a programming session, a message notifies you of the battery issue. You need to either recharge the battery in the rechargeable neurostimulator or replace the batteries in the ENS.

For information on battery service life for implanted neurostimulators, refer to "Elective Replacement Indicator (ERI) message" on page 60 and "End Of Service (EOS) message" on page 60.

Entering implant location

The implant location for implantable neurostimulators is entered on the Device screen. The Device screen is accessed through the workflow navigator.

- Select the implant location from the scrollable drop-down list.

Note: For neurostimulators with SureScan MRI Technology (97702, 97712, 97713, 97714), the implant location impacts the MRI scan eligibility shown on the MRI screen.

Entering patient information

Patient information such as patient name and diagnosis is entered on the Patient Info screen. The Patient Info screen is accessed through the workflow navigator or by tapping Patient Info on the Side Menu (≡) when the screen is not part of a workflow.

- Enter the applicable patient information. Red asterisks indicate the required fields.

Entering lead information

The Lead Select screen is used to identify which lead models are implanted, associate the lead electrode numbering to the neurostimulator electrode numbering, and check connectivity of the lead or extension contacts with the neurostimulator contacts. The Lead Select screen is accessed through the workflow navigator.

- To assign a lead, select the lead model from the scrollable drop-down list. Press the lead image, and then drag and drop it at the appropriate electrode numbers of the device.
- To delete a lead, tap the lead image, and then tap the X that appears.
- To change a lead, you must first delete the current lead, and then assign a new lead.

Notes:

- Make sure the electrode numbering on the screen matches the electrode numbering of the implanted leads. If necessary, you can swap leads using the drag-and-drop feature.
- On surgical leads, the blue box on the lead arm represents the white marker band on the actual lead. The marker band indicates the lead arm for electrodes 0-7.
- For neurostimulators with SureScan MRI Technology (97702, 97712, 97713, 97714), the assigned leads impact the MRI scan eligibility shown on the **MRI** screen.

Testing lead connectivity

Once the implanted lead(s) are assigned on the **Lead Select** screen, you can check the connectivity of the lead or extension contacts with the neurostimulator contacts. Prior to running the test, ensure the leads or extensions are properly connected to the neurostimulator.

- Tap the **CHECK CONNECTIVITY** button (🔍). (If the button is not visible, tap **UPDATE** on the **Lead Select** screen to save the lead data and activate the button.)
 - Green electrodes indicate that the lead or extension contacts have a good connection to the neurostimulator contacts.
 - Red electrodes (with an "X") indicate a bad connection.
 - If there are red electrodes, gently reconnect the leads or extensions for a better connection to the neurostimulator contacts. Then repeat the connectivity check.

Entering extension, abandoned lead, or tip location information

The **Tip Location** screen is accessed through the workflow navigator. The options on the screen will vary based on the neurostimulator model being programmed. The screen is not supported for Itrel 4 (37703, 37704).

- Check the **Yes** or **No** checkbox to indicate whether extensions are implanted.
- Check the **Yes** or **No** checkbox to indicate whether abandoned leads remain implanted (for MRI models only).
- Select the tip location for each lead from the scrollable drop-down list (for MRI models only).

Note: For neurostimulators with SureScan MRI Technology (97702, 97712, 97713, 97714), the abandoned lead, extension, and tip location information impacts the MRI scan eligibility shown on the **MRI** screen.

Clearing information in external neurostimulator for re-use

The Model 37022 External Neurostimulator is a re-usable device. The **Clear Device** screen accessed through the **Configured ENS: End Evaluation** workflow is used to clear the patient and therapy information from the device and prepare the ENS for re-use.

- On the **Clear Device** screen, tap the **Clear Device** button, and then confirm the action.
- A message appears after the device is cleared. Tap **OK** to end the session and return to the initial app screen.

The next time the external neurostimulator is interrogated, the device is identified as a new ENS and the **New ENS: Start Usage > Start Evaluation** workflow appears.

Working with Groups and Programs

Programs screen

The **Programs** screen is used to create, copy, move, and delete programs within groups and to define the active group. The **Programs** screen is accessed through the workflow navigator.

- Only one group can be active at a time. A highlighted box around a group indicates the active group.
- To make another group active, tap that group name.
- To copy or move a program, press and drag that program to the intended location, and then select **MOVE** or **COPY** from the message that appears.
- To delete a program, press and drag that program to the **Trash** icon (🗑).
- To create a new program in the active group, tap the add symbol (+) on the desired program. The **Program** subscreen will display where you can assign the program settings.
- To view the program settings for an existing program, first make the group active, and then tap the program to view the **Program** subscreen.

Note: For neurostimulators with AdaptiveStim Technology (37714, 97714), the AdaptiveStim indicator (⚡) appears next to the group name if AdaptiveStim Technology is enabled for the group.

Program subscreen

The **Program** subscreen is used to assign, review, and change program settings for programs in the active group. The **Program**

subscreen is accessed by tapping a program on the **Programs** screen.

- A line is shown under the program number that is currently displayed on the **Program** subscreen.
- To move from program to program within the active group, swipe left or right on the screen. Alternately, you can tap a program number.
- Use the **Program** subscreen to assign the electrode configuration for the leads and to set pulse width, rate, and amplitude. Refer to "Adjusting Stimulation Parameters" on page 31.
- From the **Program** subscreen, you can access the **Lead Manipulation**, **Energy**, and **Patient Access** subscreens. Refer to the following sections for details about these subscreens:
 - "Positioning Leads Relative to Implant Location" on page 34.
 - "Customizing SoftStart/Stop and Cycling" on page 35.
 - "Enabling Patient Access" on page 36.

Turning stimulation on or off

The **Stimulation** toggle is found on the action bar at the top-right corner of the screen. The toggle is available when at least one program exists. When no programs exist, the area in the action bar reads "Stimulation Not Setup".

- Tap the **Stimulation** toggle to turn stimulation on (🔘) or off (🔘) at any time.

Note: For neurostimulators with SureScan MRI Technology (97702, 97712, 97713, 97714), the MRI Status indicator (📶) displays in the action bar whenever stimulation is off. The MRI Status indicator does not provide MRI scan eligibility information.

Adjusting Stimulation Parameters

Stimulation is delivered in electrical pulses. The stimulation parameters, which include electrode configuration, pulse width, rate, and amplitude, are attributes of these pulses and can be adjusted to manage patient therapy.

High-output interlocks

This feature applies to RestoreUltra (37712, 97712), RestoreSensor (37714, 97714), and ENS (37022).

Certain combinations of high amplitude, pulse width, and rate settings are not allowed by the clinician programmer app. High-output interlocks can prevent certain values and features from being available for programming. If you attempt to program a parameter value that will cause the settings to exceed the high-output interlock limit, the desired parameter value can only be achieved by reducing one of the other parameter values.

Shared electrodes

This feature applies to RestoreUltra (37712, 97712), RestoreSensor (37714, 97714), and ENS (37022).

One electrode may be active in more than one program within a group. When more than one electrode is shared by more than one program, certain restrictions are applied to prevent rate and pulse width upper limits from exceeding 260 Hz for rate or 450 μ s for pulse width. The collective rate for a group is calculated by multiplying the rate assigned to each program by the number of programs in that group. Shared electrodes appear grayed out on the **Program** subscreen when unavailable for programming.

Assigning electrode configuration

Electrode configuration for the leads is assigned on the **Program** subscreen.

Electrode configuration must include at least one negative (–) and one positive (+) electrode. You cannot set amplitude for the program until you have assigned at least one negative (–) and one positive (+) electrode to the lead.

- To assign polarity, tap on the appropriate electrode. Tap once for negative (–) polarity, twice for positive (+) polarity, or a third time to remove polarity. Tap **Update** when finished.
- To delete an electrode, press and drag that electrode to the **Trash** icon (🗑️).
- To delete all assigned electrodes, tap the **Select All** button (⌵), and then tap the **Trash** icon (🗑️).
- To move all electrodes at once, tap the **Select All** button (⌵), and then press and drag the electrodes to a new location on a lead. This action will set the amplitude to 0 V. The move feature is not available for Itrel 4 (37703, 37704).
- To move an individual electrode, you must first delete the electrode and then assign a new electrode.

Note: Adding or deleting an electrode after the amplitude is assigned will reset the amplitude to 0 V by default.

For Itrel 4 (37703, 37704), the neurostimulator can be programmed to deliver unipolar stimulation (using the neurostimulator case as a positive electrode). To assign the case as a positive electrode, tap the 'C' on the lead image.

Assigning pulse width, rate, and amplitude

Pulse width, rate, and amplitude are assigned on the **Program** subscreen.

During a programming session, the selected values for pulse width, rate, and amplitude are automatically programmed into the neurostimulator. Once amplitude is set and stimulation is turned on, the neurostimulator changes stimulation to the new setting every time a stimulation parameter is changed.

△ **Caution:** To prevent possible uncomfortable or unexpected stimulation (jolting or shocking sensation) during a programming session, adjust stimulation parameters in small increments above the perception threshold (the parameter values at which the patient first perceives a sensation of stimulation).

△ **Caution:** Changing the pulse width or rate in a program may affect AdaptiveStim therapy. Test each of the AdaptiveStim postures in the group for appropriate stimulation after making pulse width or rate changes. Refer to "Testing AdaptiveStim settings" on page 48.

Assigning pulse width

Pulse width is the duration of the pulses in microseconds (µs).

△ **Caution:** To prevent possible uncomfortable or unexpected stimulation (jolting or shocking sensation), decrease the amplitude to the perception threshold (the amplitude at which the patient first perceives a sensation of stimulation) before changing the pulse width. After changing the pulse width, slowly increase the amplitude.

1. In the **Program** subscreen, tap **PULSE WIDTH** to access the pulse width control.
2. Select the pulse width using any of these controls:
 - Scrollable drop-down list.

- Dial (press and drag the dot on the dial).
 - Back arrows (◀) to decrease and forward arrows (▶) to increase.
3. If stimulation is on, the **Stop** and **Jump** buttons become available when the pulse width is increased. Use these buttons, as needed.
 - Tap the **Stop** button (□) to stop stimulation from increasing to the new value.
 - Tap the **Jump** button (▶▶) to jump to the new value without incremental increases.

Assigning rate

Rate is the frequency of pulses per second, in hertz (Hz). For the devices that support more than one program in a group, the group rate is the rate for each program in that group.

1. In the **Program** subscreen, tap **GROUP RATE** to access the rate control.
2. Select the rate using any of these controls:
 - Scrollable drop-down list.
 - Dial (press and drag the dot on the dial).
 - Back arrows (◀) to decrease and forward arrows (▶) to increase.

Assigning amplitude

Amplitude is the strength of the pulse in volts (V). The electrode configuration for the leads must be assigned before you can set amplitude.

1. In the **Program** subscreen, tap **AMPLITUDE** to access the amplitude control.

2. Select the amplitude using any of these controls:

- Scrollable drop-down list.
- Dial (press and drag the dot on the dial).
- Back arrows (◀) to decrease and forward arrows (▶) to increase.

Note: Amplitude increases in increments. When amplitude is decreased, amplitude decreases immediately to the new value.

3. Use these buttons, as needed:

- Tap the **Stop** button (□) to stop stimulation from increasing to the new value.
- Tap the **Jump** button (▶▶) to jump to the new value without incremental increases.
- Tap the **Amplitude to zero** button (⊖) to immediately set amplitude to 0 V.

Positioning Leads Relative to Implant Location

Use the **Lead Manipulation** subscreen (accessed from the **Program** subscreen) to represent the relative positioning of the implanted leads as viewed from the patient's back. The subscreen is only applicable if more than one lead is configured.

- To swap the relative position of the leads (for example, to switch the position of lead #2 with lead #1), press and drag a lead to the location of the other lead.

The positioning of the leads replaces how the leads are shown on the **Programs** screen and **Program** subscreen.

Note: The **Lead Manipulation** screen is not supported for Itrel 4 (37703, 37704) or for surgical leads.

Customizing SoftStart/Stop and Cycling

Use the **Energy** subscreen (accessed from the **Program** subscreen) to customize the optional device settings for the SoftStart/Stop and Cycling features. These features are designed to increase patient comfort and ease of use.

Changing SoftStart/Stop settings

SoftStart/Stop slowly increases the amplitude when stimulation is turned on and slowly decreases the amplitude when stimulation is turned off. The slow ramping may feel more comfortable to sensitive patients.

SoftStart/Stop is assigned at the device level: that is, it is either on or off for all groups. The SoftStart/Stop feature defaults to "on" with a duration of 4 seconds.

- In the **Energy** subscreen, use the **SOFTSTART/STOP** drop-down list to select a duration or to turn SoftStart/Stop off.

If AdaptiveStim Technology and SoftStart/Stop are both enabled and the amplitude increases or decreases because of a change in position:

- The amplitude will slowly increase using the ramping programmed for SoftStart/Stop.
- The amplitude will decrease immediately without ramping.

Turning cycling on or off

Cycling turns stimulation on and off at clinician-determined intervals. Due to a carryover effect, the patient may continue to experience symptom suppression during the cycling off time. The default setting for the cycling feature is "off".

Refer to the *System Eligibility, Battery Longevity, Specifications reference manual* for information about how cycling affects battery longevity and recharge intervals.

1. In the **Energy** subscreen, tap the **CYCLING** toggle to turn cycling on () or off (.
2. If cycling is on, select the durations for the **ON TIME** and **OFF TIME** intervals.
3. Tap the group name toggle to turn cycling on () or off () for a specific group. The setting affects all programs in the group.

Note: Cycling is not available when AdaptiveStim Technology is enabled for a group.

Enabling Patient Access

Use the **Patient Access** subscreen (accessed from the **Program** subscreen) to assign patient controls, if appropriate for the patient. There are three patient controls on the **Patient Access** subscreen: **Limits**, **TargetMyStim**, and **Adjustment**.

Note: The patient controls may vary depending on neurostimulator model.

Limits

The **Limits** patient control on the **Patient Access** subscreen provides the ability to enable or disable patient control of stimulation parameters (amplitude, pulse width, group rate). When patient control is enabled for a stimulation parameter, the patient can adjust that stimulation parameter using the patient programmer.

If patient control is enabled, the upper patient limit for the stimulation parameter must be assigned. The lower patient limit is set to the lowest possible value by default and is not customizable. Patient control of stimulation parameters is assigned at the group level.

Note: By default, patient control of amplitude is enabled, and patient control of pulse width and group rate is disabled.

1. Next to the **LIMITS** heading, make sure that the group you intend to configure is shown.
2. Tap the toggle next to the stimulation parameter name to enable (☑) or disable (☐) patient control of the stimulation parameter.
3. For enabled stimulation parameters, set the upper patient limit using any of these controls:
 - Scrollable drop-down list.
 - Dial (press and drag the dot on the dial).

4. After the patient controls are defined, tap **UPDATE**.
5. To see the patient control settings for all groups, tap the **OVERVIEW** button (☰).

Amplitude

- The lowest allowable option for the upper patient limit is determined by the highest amplitude setting in the group's existing programs.
- When AdaptiveStim Technology is enabled for a group, patient control is automatically enabled for amplitude and cannot be disabled for that group.

Pulse width

- The lowest allowable option for the upper patient limit is determined by the highest pulse width setting in the group's existing programs.
- When AdaptiveStim Technology is enabled for a group, patient control is not available for pulse width.

Group rate

- The lowest allowable option for the upper patient limit is determined by the current rate setting for the group.
- When AdaptiveStim Technology is enabled for a group, patient control is not available for group rate.

Effect of subsequent parameter changes

- During a subsequent programming session, if a parameter is increased above the upper patient limit, the upper patient limit automatically changes to that new value.

TargetMyStim

Enabling the TargetMyStim patient control on the **Patient Access** subscreen allows the patient to move stimulation settings from electrode to electrode up and down the lead.

1. Under the **TARGETMYSTIM** heading, make sure that the group you intend to configure is shown.
2. Tap the group name toggle to enable  or disable  the patient control.

Notes:

- The TargetMyStim patient control is not supported for Itrel 4 (37703, 37704), Restore (37711), and RestorePrime (37701).
- TargetMyStim is only available if the following criteria are met:
 - The electrode pattern for at least one program in the group can be moved up or down the lead.
 - Amplitude limit is enabled.
 - For RestoreSensor (37714, 97714), RestoreUltra (37712, 97712), and ENS (37022): No program in the group has a pulse width greater than 450 µs; the collective rate of all programs in the group is not greater than 260 Hz; there are no shared electrodes.
- The TargetMyStim patient control is not available when AdaptiveStim Technology is enabled for a group.

Adjustment

Enabling the Adjustment patient control on the **Patient Access** subscreen allows the patient to increase or decrease all program amplitudes in the active group at the same time.

- Under the **ADJUSTMENT** heading, tap the **ALL PROGRAMS TOGETHER** toggle to enable  or disable  the patient control.

Note: The Adjustment patient control is not supported for Itrel 4 (37703, 37704), Restore (37711), and RestorePrime (37701).

Programming Scheduled Therapy

The Scheduled Therapy feature is used to activate specific groups at programmed times during a 24-hour period. This feature is not supported for Itrel 4 (37703, 37704) and RestoreSensor (37714, 97714).

Therapy groups can be scheduled to activate and deactivate at specific times during a 24-hour period. Each activation and deactivation of a group is called an event. Up to 8 scheduled events can be programmed during a 24-hour period. Events cannot overlap, and only one group can be active at any given time. The minimum time for an event is 30 minutes.

The **Schedule** screen is accessed from the Side Menu (). The groups available for scheduling are shown in the **GROUP** list. When events are scheduled, they display in the **DAILY SCHEDULE**.

- To schedule a new event, tap on a group in the **GROUP** list and then assign the start and end times in the pop-up window.
(Alternately, you can drag a group from the **GROUP** list and drop it on the **DAILY SCHEDULE**. Then adjust the start and end times in the pop-up window.)
- To edit an existing event, tap the event in the **DAILY SCHEDULE**, and then adjust the start and end times in the pop-up window.
(Alternately, you can move an event in the **DAILY SCHEDULE** using the drag-and-drop feature.)
- To remove an existing event, drag the event from the **DAILY SCHEDULE** and drop it on the **Trash** icon ().
(Alternately, you can tap an event in the **DAILY SCHEDULE**, and then tap the **Trash** icon ()) in the pop-up window.)

- To turn the Scheduled Therapy feature on or off, use the on/off toggle (green indicates that Scheduled Therapy is on).

Note: Scheduled Therapy operation will be suspended during a programming session to prevent therapy changes from occurring during the session. When the session is exited, Scheduled Therapy operation resumes.

Checking System Performance

The measurement functions on the **Electrode Impedance** screen and **Group Impedance** screen assist in identifying problems with system components or the entire implanted system.

Measurements and diagnostic data obtained from the clinician programmer app are intended to aid in your clinical assessment. However, as with any electronic system, internal and external factors can influence neurostimulator measurements. For example, changes in lead position can affect the stimulation current or the impedance measurement. If you obtain a reading that seems inconsistent with your observations, repeat the measurement.

Before using any measurement functions, please note the following:

- **Measure electrode impedance at the beginning of each programming session.** These measurements verify the integrity of lead, extension, and connector pathways. They may provide information about lead problems (such as lead breakage, short circuit, open circuit). For example, measurements that show a significant increase in electrode impedance can indicate a fractured lead conductor or a loose setscrew. Conversely, a significant decrease in electrode impedance can indicate shorted conductors or a break in lead insulation. Measurements taken at the beginning of the session may be useful in interpreting diagnostic data collected since the previous session.
- **Measure group (therapy) impedance at the end of each programming session.** These measurements provide documentation that pathways are intact and the current provided is sufficient for the selected therapy.

Testing electrode impedance

The **Electrode Impedance** screen is accessed through the workflow navigator or by tapping **Impedance** on the Side Menu (≡) when not part of a workflow.

For Restore (37711) and RestorePrime (37701), preset rate, pulse width, and amplitude values are used during the electrode impedance test. For all other models, preset rate and pulse width values are used, and the amplitude value is selectable on the **Start Measurement** screen. At the end of the test, the most recently programmed settings are restored.

1. On the **Electrode Impedance** screen, tap the **MEASURE ALL IMPEDANCES** button (⚡).
2. On the **Start Measurement** screen, select the amplitude value, if applicable.
3. Tap the **START** button.
4. After the test is complete, review the impedance results.

Notes:

- If all impedance values are within range, the message "All impedances are within normal range" appears in the **GOOD** box.
- If any impedance value is out of range (low or high), the results appear in the **OUT OF RANGE** box. Refer to "Troubleshooting electrode impedance results" on page 40.
- The impedance value shown for each electrode is measured with respect to the selected reference electrode. Use the **REFERENCE ELECTRODE** scrollable drop-down list to select a different reference electrode, if applicable.
- Once an impedance test is performed, results can be viewed at any time during the programming session by returning to the

Electrode Impedance screen. You can also view or print the Session Report, which shows the impedance measurements.

Troubleshooting electrode impedance results

The OUT OF RANGE box on the Electrode Impedance screen identifies where low impedance or high impedance was measured in the system. The example in Figure 4 shows both low impedance and high impedance results.

< 50 Ω
13–15
> 10000 Ω
0: 14
1: 14
2: 14
3: 14
4: 14
5: 14
6: 14
7: 14
8: 14
9: 14
10: 14
11: 14
12: 14
13: 14
14: 0, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 15
15: 14

Figure 4. Example of out-of-range electrode impedance results.

In this example, a low impedance value was measured for the electrode pair 13–15, and high impedance values were measured for several electrode pairs with electrode 14 common to those pairs. Electrode 14 is likely the electrode causing the high impedance issue.

Refer to the troubleshooting details for high impedance and low impedance on the following pages. For further assistance with

troubleshooting electrode impedance results, contact Medtronic using the contact information listed on the back cover of this manual.

High impedance troubleshooting

High impedance suggests a possible open circuit in the system.

High impedance for Restore (37711) and RestorePrime (37701):
>3,600 ohms

High impedance for all other models:

- At 0.25 V: >4,000 ohms
- At 0.70 V: >10,000 ohms
- At 1.50 V: >20,000 ohms
- At 3.00 V: >40,000 ohms

Follow these steps to troubleshoot high impedance values:

1. If applicable, perform the test again at the next highest amplitude up to 3.00 V (if within the patient's comfort range).
2. If impedance values remain high, consider these troubleshooting options:
 - Assess whether therapy can be achieved using other electrodes.
 - Check for misprogrammed information. For example, check the **Lead Select** screen to ensure the electrode numbering on the screen matches the electrode numbering of the implanted lead(s).
 - Consider using x-ray to inspect the system components for disconnects or damage.
3. If troubleshooting measures fail to resolve the high impedance values, consider replacing system components as needed.

Low impedance troubleshooting

Low impedance suggests a possible short circuit in the system that can lead to rapid battery depletion or inadequate therapy. Do not use electrodes with a low impedance.

Low impedance for Restore (37711) and RestorePrime (37701):
<70 ohms

Low impedance for all other models:

- At 0.25 V: <50 ohms
- At 0.70 V: <50 ohms
- At 1.50 V: <50 ohms
- At 3.00 V: <150 ohms

Testing group impedance

The **Group Impedance** screen is accessed through the workflow navigator or by tapping **Impedance** on the Side Menu (≡) when not part of a workflow.

Group (therapy) measurements are impedance and stimulation current measurements taken at the programmed settings for each program in a selected group.

 **Caution:** During the group impedance test, stimulation changes to the parameter values for each program as the impedance is measured for the group. This may affect patient stimulation while the test is in progress. You may need to decrease stimulation parameters to avoid patient discomfort.

1. On the **Group Impedance** screen, tap the group name you want to test, and then tap the **MEASURE IMPEDANCES** button (Ω).

2. The stimulation parameter ranges to be tested display on the **Start Measurement** screen. Tap the **START** button after confirming the values.
3. Review the impedance and stimulation current results. Refer to "Troubleshooting group impedance results" on page 41.

Notes:

- The neurostimulator reports current measurement values as current delivered during the stimulation pulse.
- Once an impedance test is performed, results can be viewed at any time during the programming session by returning to the **Group Impedance** screen. You can also view or print the Session Report, which shows the impedance measurements.

Troubleshooting group impedance results

When impedance results are in the normal range, the impedance value (in ohms) and the stimulation current (in milliamps) display on the **Group Impedance** screen. Normal impedance ranges are:

- 70 to 3,600 ohms for Restore (37711) and RestorePrime (37701)
- 50 to 4,000 ohms for all other models

When impedance results are outside of the normal range, an indicator displays (either < ## Ω or > #### Ω) where the value is the low or high range threshold for the neurostimulator. Perform an electrode impedance test to help determine the likely cause of the circuit problem. Refer to "Testing electrode impedance" on page 39.

Notes:

- If an impedance result is not available or out-of-regulation occurs, the message "HIGH AMPLITUDE" displays. Decrease

stimulation parameters or reduce the number of electrodes for that program, and then repeat the measurement.

- If a program amplitude is too low for testing, the message "LOW AMPLITUDE" displays. Increase the amplitude for that program, and then repeat the measurement. Minimum amplitude values are:
 - 0.24 V for Restore (37711) and RestorePrime (37701)
 - 0.25 V for all other models

Preparing System for an MRI Scan

Refer to the *MRI Guidelines for Medtronic Neurostimulation Systems for Chronic Pain instructions for use* manual for the MRI conditions and MRI-specific warnings and cautions for conducting an MRI scan.

- Stimulation must be off during an MRI scan.
- To turn stimulation off, use the **Stimulation** toggle () on the action bar. Gray indicates that stimulation is off.

View MRI workflow

The View MRI workflow is only applicable to neurostimulators with SureScan MRI Technology (97702, 97712, 97713, 97714).

1. On the **Device**, **Lead Select**, and **Tip Location** screens, confirm that the following component information is complete and accurate:
 - Neurostimulator and lead model numbers
 - Neurostimulator implant location and lead tip location
 - Presence of extensions or abandoned leadsThis information impacts the MRI scan eligibility shown on the **MRI** screen.
2. On the **Patient Info** screen, confirm that the patient information is complete and accurate.
3. On the **MRI** screen, review the MRI scan eligibility.
4. Determine whether the patient is eligible for the intended MRI scan. Do not proceed if eligibility cannot be determined. Contact your Medtronic representative, if needed.

5. To prepare the device for an MRI scan, turn stimulation off using the **Stimulation** toggle () on the action bar. Gray indicates that stimulation is off.
6. Confirm that the message "Device ready for MRI scan" displays in the status area and the MRI Status indicator () displays in the action bar.
7. On the **Reports** screen, select **MRI Report** from the drop-down list. Refer to "Working with Reports" on page 55 for details about viewing and printing reports.

Note: If printing capabilities are not available, you can manually complete the MRI Eligibility Form in the *MRI Guidelines for Medtronic Neurostimulation Systems for Chronic Pain instructions for use* manual.

8. Inform the patient of the following:
 - Stimulation must be off during the MRI scan.
 - Because stimulation is off, pain symptoms may return.
 - Bring the most up-to-date patient ID card to all MRI appointments.
 - Bring the patient programmer to all MRI appointments.
 - After the MRI scan is complete and the patient is outside of the magnet room, stimulation can be turned back on using the patient programmer. Or, if appropriate, the patient can revisit the clinician after the MRI scan to have stimulation turned back on with the clinician programmer app.
9. Give the MRI Report to the patient or send the report to the MRI center.

The MRI clinician uses the MRI Report to confirm the neurostimulator model, session date, patient information, and MRI scan eligibility of the patient's implanted system.

AdaptiveStim Technology

The RestoreSensor neurostimulators with AdaptiveStim Technology (37714, 97714) can be programmed to adjust amplitude automatically when the patient changes body position. Amplitude settings can be assigned to the following positions:

- Upright
- Mobile
- Lying Front
- Lying Back
- Lying Right
- Lying Left

When AdaptiveStim Technology is turned on for a patient, the neurostimulator automatically adjusts the amplitude based on the patient's body position and the clinician-programmed settings. If the patient makes an amplitude adjustment with the patient programmer, the neurostimulator associates the new amplitude with the patient's body position as long as the patient stays in that position for at least three minutes. The next time the patient enters that position, the neurostimulator will return to the new amplitude. In this way, AdaptiveStim Technology can respond to patient feedback.

To accommodate the individual needs of the patient, the following AdaptiveStim settings can be adjusted by the clinician:

- **Position angles:** The range or area of a position.
- **Transition times:** The amount of time that must elapse after moving into a new position (and remaining in the new position) before the amplitude will change.
- **Mobility rate:** The intensity of movement that is considered a change from the stationary position of Upright to the moving position of Mobile.

Caution regarding cervical location

△ **Caution:** The use of AdaptiveStim Technology associated with position changes in the cervical location has not been shown to be safe and effective. The use of AdaptiveStim Technology associated with position changes for locations outside of the thoracic or lumbar spine may result in unintended stimulation.

Caution regarding healing period

△ **Caution:** Do not configure AdaptiveStim Technology during the implant procedure. A four-week healing period is recommended after implanting the neurostimulator to allow time for the implanted components to stabilize. If the tissue around the neurostimulator is not sufficiently healed, movement of the neurostimulator is possible and may affect orientation.

Configuring AdaptiveStim Technology with the wizard

When you configure AdaptiveStim Technology with the setup wizard, all programmed groups and every position including Mobile are AdaptiveStim-enabled, no matter what settings were previously assigned to those groups.

When a group is AdaptiveStim-enabled, the AdaptiveStim indicator (i) appears next to the group name on the **Programs** screen, and the following programming features are not available for that group:

- Patient access for rate
- Patient access for pulse width
- Cycling
- TargetMyStim

Note: Amplitude limit is automatically enabled when AdaptiveStim Technology is configured for a group and cannot be disabled for that group.

1. Go to the **Programs** screen, and tap a group name to activate the group you want to configure.
2. Select **AdaptiveStim** in the workflow, and tap **CONFIGURE** to start the setup wizard.
3. Follow the on-screen instructions for all applicable positions.
 - On the orientation screen, have the patient get in the position and tap **NEXT** to complete the orientation.
Note: You can tap **SKIP** to bypass certain positions. Orientation still occurs for the position and a position amplitude of zero is assigned to that position.
 - On the programming screen, tap the program number you want to configure, and then set the position amplitude using the amplitude control. Refer to "Assigning amplitude" on page 32.
 - If applicable, tap another program number, and then set the position amplitude. When finished, tap **NEXT** to go to the next position.
4. When the wizard is finished, a confirmation screen displays stating that you have successfully configured AdaptiveStim Technology for the group.
5. Tap **CLOSE** to go to the **AdaptiveStim** screen.

Reconfiguring an AdaptiveStim-enabled group

If you want to reconfigure an AdaptiveStim-enabled group, complete these steps to delete the existing AdaptiveStim settings and reassign

position amplitudes for that group. You also have the option to reorient patient positions.

1. Go to the **Programs** screen, and tap a group name to activate the group you want to reconfigure.
2. Go to the **AdaptiveStim** screen, tap the **Settings** button (⚙️), and then tap **Configure**.
3. Tap **Start Over**.
4. (Optional) Tap the **Reorient** checkbox if you also want to reorient the neurostimulator to the patient positions.
5. Tap **Yes** to proceed.
6. Follow the on-screen instructions to reorient the neurostimulator to the patient positions (if applicable) and reassign position amplitudes for the programs in the active group.
7. When the wizard is finished, a confirmation screen displays stating that you have successfully reconfigured AdaptiveStim Technology for the group.
8. Tap **CLOSE** to go to the **AdaptiveStim** screen.

What you can do on the AdaptiveStim screen

The **AdaptiveStim** screen displays when an AdaptiveStim-enabled group is the active group and you select **AdaptiveStim** in the workflow. The **AdaptiveStim** screen can also be accessed by tapping **AdaptiveStim** on the Side Menu (≡) when not part of a workflow. Refer to Table 11 on page 46 for information about the tasks you can complete on the **AdaptiveStim** screen.

Table 11. Tasks on the AdaptiveStim screen.

Tasks on the AdaptiveStim screen	More information
Review position amplitudes for all AdaptiveStim-enabled groups. Use this task to view a summary of the position amplitudes for all AdaptiveStim-enabled groups and programs.	Page 48
Test AdaptiveStim settings. Use this task to test AdaptiveStim settings during a programming session.	Page 48
Adjust position amplitudes for individual positions. Use this task to adjust position amplitudes for individual positions in the active group.	Page 48
Check patient positions. Use this task to determine if the neurostimulator correctly detects the patient positions.	Page 49
Reorient patient positions. Use this task to reorient the neurostimulator to the patient positions without changing any group settings.	Page 49
Reconfigure an AdaptiveStim-enabled group. Use this task to delete all existing AdaptiveStim settings for a group and reassign position amplitudes. You also have the option to reorient patient positions.	Page 45
Disable AdaptiveStim Technology for a group. Use this task to clear all existing AdaptiveStim settings for a group. All other AdaptiveStim-enabled groups are not impacted.	Page 49
Enable or disable AdaptiveStim Technology for the patient. If enabled, AdaptiveStim Technology and all assigned settings will be on when the programming session ends. If disabled, AdaptiveStim Technology and all assigned settings will be off when the programming session ends. However, the patient can turn on AdaptiveStim Technology using the patient programmer.	Page 49
Turn the Position Trend diary on or off. The Position Trend diary shows the time the patient spent in each position during the periods between the last four programming sessions.	Page 50
Adjust position angles. Position angle is the range or area of a position and is used to ensure the neurostimulator properly recognizes the Upright position from the Lying positions.	Page 50
Adjust transition times. Transition time is the amount of time that must elapse between changes in position before the amplitude will change.	Page 50

Tasks on the AdaptiveStim screen	More information
Adjust mobility rate. Mobility rate is the intensity of movement that is considered a change from the stationary position of Upright to the moving position of Mobile.	Page 51
Turn off mobility detection. Turning off mobility detection disables the Mobile position for AdaptiveStim Technology.	Page 51
Troubleshoot AdaptiveStim Technology. Table 12 lists potential issues that can occur with AdaptiveStim Technology, as well as possible solutions for resolving those issues.	Page 51

Reviewing position amplitudes for all AdaptiveStim-enabled groups

A summary of the position amplitudes for all AdaptiveStim-enabled groups and programs can be viewed from the **AdaptiveStim** screen.

1. Go to the **AdaptiveStim** screen, tap the **Settings** button (⚙️), and then tap **Overview**.
2. Review the position amplitudes assigned to the AdaptiveStim-enabled groups. If all groups and programs are not visible, you may need to swipe left or right to see all of the data.
3. Tap **CLOSE** to return to the **AdaptiveStim** screen.

Testing AdaptiveStim settings

During a programming session, AdaptiveStim settings are suspended. The Test AdaptiveStim feature can be used to test AdaptiveStim settings in session and obtain patient feedback. There must be at least two different position amplitudes defined in a program to run the test.

1. Go to the **Programs** screen, and tap a group name to activate the group you want to test.
2. Go to the **AdaptiveStim** screen, and tap the **TEST ADAPTIVESTIM** button (🧪).
3. Follow the on-screen instructions to test the AdaptiveStim settings and obtain patient feedback.
4. Tap **STOP** to end the test.

Adjusting position amplitudes for individual positions

Complete these steps to adjust position amplitudes for individual positions in a group.

1. Go to the **Programs** screen, and tap a group name to activate the group you want to adjust.
2. Go to the **AdaptiveStim** screen, and tap the position you want to adjust.

On the screen, you see the current programs in the group. Within each program box, the program amplitude is displayed in white and the position amplitude is displayed in blue. A white box around a program identifies the active program.

3. Tap the program number you want to adjust.
4. Use the amplitude control to adjust the program amplitude (the stimulation the patient is receiving). Refer to "Assigning amplitude" on page 32.

The white dot represents the program amplitude. As you make amplitude changes, the white dot moves to reflect the current amplitude.

Note: You can also drag the white dot to assign amplitude. Stimulation ramps up to the amplitude value.

5. When you find an acceptable amplitude, tap **SET POSITION AMPLITUDE** to assign the stimulation the patient is receiving as the position amplitude.
The blue dot represents the position amplitude.
6. When finished making adjustments, tap **CLOSE** to return to the **AdaptiveStim** screen.

Checking patient positions

Complete these steps to determine if the neurostimulator correctly detects the patient positions.

1. Place the patient in the position you intend to check.
2. Go to the **AdaptiveStim** screen, and tap the **CHECK POSITION** button (🔄).
3. Confirm that the position shown on the **AdaptiveStim** screen corresponds to the patient's body position.

Note: If the patient is in a position between the Upright and Lying positions, the app identifies the position as Transition Zone.

Reorienting patient positions

Complete these steps to reorient the neurostimulator to the patient positions without changing any group settings.

1. Go to the **AdaptiveStim** screen, tap the **Settings** button (⚙️), and then tap **Configure**.
2. Tap **Reorient Only**.
3. Tap **Yes** to clear the current orientation.
4. Follow the on-screen instructions to reorient the neurostimulator to the patient positions.
5. (Optional) Tap **SKIP** if you want to bypass a position.
6. When all positions have been reoriented, a confirmation screen displays.
7. Tap **CLOSE** to go to the **AdaptiveStim** screen.

Disabling AdaptiveStim Technology for a group

Complete these steps to clear all existing AdaptiveStim settings for a group. All other AdaptiveStim-enabled groups will not be impacted.

1. Go to the **Programs** screen, and tap a group name to activate the group.
2. Go to the **AdaptiveStim** screen, and tap the on/off toggle next to the **Settings** button (⚙️).
 - Green indicates that the group is enabled.
 - Gray indicates that the group is disabled.
3. Tap **UPDATE**, and then confirm the disable notification.

Enabling or disabling AdaptiveStim Technology for the patient

You can choose to either enable or disable AdaptiveStim Technology for the patient. If enabled, AdaptiveStim Technology and all assigned settings will be on when the programming session ends. If disabled, AdaptiveStim Technology and all assigned settings will be off when the programming session ends. However, the patient can turn on AdaptiveStim Technology using the patient programmer.

1. Go to the **AdaptiveStim** screen, and tap the **RESUME** on/off toggle.
 - Green indicates that AdaptiveStim is enabled.
 - Gray indicates that AdaptiveStim is disabled.
2. Tap **UPDATE** to accept the changes.

Turning the Position Trend diary on or off

The Position Trend diary shows the time the patient spent in each position during the periods between the last four programming sessions. The diary is accessed through the **Diaries** screen.

By default, the Position Trend diary is turned on when AdaptiveStim Technology is initially configured or reconfigured, or when patient positions are reoriented. If the diary is turned off, from that point, position data is not saved. If the diary is on and AdaptiveStim Technology is turned off, position data continues to be tracked.

1. Go to the **AdaptiveStim** screen, tap the **Settings** button (⚙️), and then tap **Settings**.
2. On the **Options** tab, tap the on/off toggle. Green indicates that the diary is on.
3. Tap **UPDATE** to accept the changes.
4. Tap **CLOSE** to return to the **AdaptiveStim** screen.

Adjusting position angles

Position angle is the range of a position and is used to ensure the neurostimulator properly recognizes the Upright position from the Lying positions. By default, the Upright position angle is set to X-Small and the Lying position angle is set to Large. All lying postures (Back, Front, Left, Right) use the same position angle value. There is a Transition Zone between the Upright and Lying positions. The size of the Transition Zone is calculated based on the Upright and Lying position angles.

1. Go to the **AdaptiveStim** screen, tap the **Settings** button (⚙️), and then tap **Settings**.
2. Tap the **Angles** tab.

3. On the **Angles** diagram, tap the position name for the angle you want to adjust (Upright, Lying Back, or Lying Front).
4. Press and drag one of the delimiter dots to set a larger or smaller angle.
5. Make other adjustments to the position angles as needed.
6. Tap **UPDATE** to accept the changes.
7. Tap **CLOSE** to return to the **AdaptiveStim** screen.

Note: To test the changes, refer to "Testing AdaptiveStim settings" on page 48.

Adjusting transition times

Transition time is the amount of time that must elapse between changes in position before the amplitude will change. You can adjust transition times for each position change.

1. Go to the **AdaptiveStim** screen, tap the **Settings** button (⚙️), and then tap **Settings**.
2. Tap the **Transitions** tab.
3. Swipe up and down on the screen to scroll and find the position change you want to adjust.
4. Tap the scrollable drop-down list and select a duration.
5. Make adjustments to other transition times as needed.
6. Tap **UPDATE** to accept the changes.
7. Tap **CLOSE** to return to the **AdaptiveStim** screen.

Note: To test the changes, refer to "Testing AdaptiveStim settings" on page 48.

Adjusting mobility rate

Mobility rate is the intensity of movement that is considered a change from the stationary position of Upright to the moving position of Mobile. The mobility rate is only accessible if the Mobile posture is enabled.

1. Go to the **AdaptiveStim** screen, tap the **Settings** button (⚙️), and then tap **Settings**.
2. Tap the **Mobility Rate** tab.
3. Confirm that the on/off toggle is on. Green indicates that the Mobile posture is enabled.
4. Press and drag the **Mobility Rate** dot on the slider.
 - **Low (L)**: Takes less motion for the neurostimulator to detect the mobile position.
 - **High (H)**: Takes more motion for the neurostimulator to detect the mobile position.
5. Tap **UPDATE** to accept the changes.
6. Tap **CLOSE** to return to the **AdaptiveStim** screen.

Note: To test the change in mobility rate, refer to "Testing AdaptiveStim settings" on page 48. To check patient positions, refer to "Checking patient positions" on page 49.

Turning off mobility detection

Turning off mobility detection disables the Mobile position for AdaptiveStim Technology. When off, the Mobile position button is dimmed in the **AdaptiveStim** screen.

1. Go to the **AdaptiveStim** screen, tap the **Settings** button (⚙️), and then tap **Settings**.

2. Tap the **Mobility Rate** tab.
3. Tap the on/off toggle. Gray indicates that the Mobile position is disabled.
4. Tap **UPDATE** to accept the changes.
5. Tap **CLOSE** to return to the **AdaptiveStim** screen.

Troubleshooting AdaptiveStim Technology

Refer to Table 12 on page 52 for a list of potential issues that can occur with AdaptiveStim Technology, as well as possible solutions for resolving those issues.

Table 12. Troubleshooting AdaptiveStim Technology.

Problem	Possible Solutions
The delay between a position change and the amplitude change is too slow or too fast for the patient.	■ Use the Transitions setting to increase or decrease the delay between position changes and amplitude changes. ■ Use the Test AdaptiveStim feature to obtain patient feedback and make adjustments as needed.
The amplitude changes too often for the patient.	■ Use the Angles setting to increase the angle for the position that is changing too often. - Or ■ Use the Transitions setting to increase the delay between position change and amplitude change.
The amplitude changes to the amplitude for Mobile while the patient is in the Upright position but not mobile.	■ Use the Mobility Rate setting to increase the intensity of movement that is required to change to the amplitude for the Mobile position. ■ Have the patient remain still in the Upright position for at least one minute. ■ Use the Check Position feature to confirm the patient's position. If "Mobile" is shown instead of "Upright" for the patient's current position, increase the mobility rate.
The amplitude does not change to the amplitude for Mobile when the patient is moving.	■ Use the Mobility Rate setting to decrease the intensity of movement that is required to change to the amplitude for the Mobile position. ■ Have the patient walk down the hall for at least one minute. ■ Use the Check Position feature to confirm the patient's position. If "Upright" is shown instead of "Mobile" for the patient's current position, decrease the mobility rate.
The patient does not like that the amplitude changes to the Mobile amplitude when starting to walk from the Upright position.	■ Turn off the Mobility Rate feature. ■ When the mobility rate is turned off, the Mobile position is no longer enabled for AdaptiveStim Technology.
Neurostimulator may have rotated or flipped within the pocket.	■ Reorient the neurostimulator to the patient positions.
Patient feels the amplitude for Upright when lying down.	■ Use the Angles setting to increase the angle for the Lying positions.
Patient is not getting adequate therapy in a position programmed with AdaptiveStim.	■ Review position amplitudes, and check that the position amplitude is not set to 0 V.

Problem	Possible Solutions
Patient unable to adjust their rate.	■ Create multiple identical groups by copying the current AdaptiveStim-enabled group. Adjust the rate for each group to provide the patient with different rates with AdaptiveStim enabled. Provide the patient guidance on how to change groups.

Viewing Diaries

The charts accessed from the **Diaries** screen provide overviews of how the patient is using the neurostimulator. The **Diaries** screen is not supported for Itrel 4 (37703, 37704).

Stimulation Usage

The **Stimulation Usage** chart shows the patient's stimulation usage for the last 30 days. The chart displays the stimulation on-time percentage for each day. The metrics provide the average percentage of stimulation usage over the last 30 days and the percentage of stimulation on-time since the last programming session.

Note: If 30 days of stimulation usage data is not available, the calculations are based on the available number of days.

Group Usage

The **Group Usage** chart displays the percentage of time that each group was used for the last 30 days or since the last programming session if less than 30 days.

Position Trend

The **Position Trend** chart is only applicable to neurostimulators with AdaptiveStim Technology (37714, 97714). The chart shows the time the patient spent in each position during the periods between the last four programming sessions. The default view shows all positions and the amount of time in percentages.

- Select a position from the scrollable drop-down list if you want to view a specific position only.
- Select **Time** if you want to see the average amount of time per day in hours and minutes.

Adjustments

The **Adjustments** chart is only applicable to neurostimulators with AdaptiveStim Technology (37714, 97714). The chart shows the number of amplitude adjustments made by the patient since the last programming session. Data is shown for each position and AdaptiveStim-enabled group.

Resting Trend

The **Resting Trend** chart is only applicable to neurostimulators with AdaptiveStim Technology (37714, 97714). The chart shows the average number of position changes from one lying position to another during the periods between the last four programming sessions.

Recharge

The **Recharge** chart is only applicable to rechargeable neurostimulators (37711, 37712, 37713, 37714, 97712, 97713, 97714). The chart displays the following information for the last six recharge sessions:

- The percentage of battery charge at the start of the recharge session
- The percentage of battery charge at the end of the recharge session
- The length of the recharge session
- The recharge session date

The metrics provide the average ending percentage of battery recharge, the average battery recharge time, the average battery recharge interval, and the median battery recharge coupling value.

Working with Reports

The following report types are available from the **Reports** screen:

- **Session Report:** Contains information about the settings programmed into the neurostimulator during a programming session as well as therapy usage and system status information.
- **MRI Report:** Applicable to neurostimulators with SureScan MRI Technology (97702, 97712, 97713, 97714). Contains information intended for MRI clinicians about the MRI scan eligibility of the implanted system for the selected patient.
- **Medtronic Data Report:** Contains information about the Medtronic data collected from the neurostimulator. This report is used by Medtronic personnel when troubleshooting. To compile a Medtronic Data Report, refer to "Compiling a Medtronic Data Report" on page 58.

The **Reports** screen is accessed through the workflow navigator or by tapping **Reports** on the Side Menu (≡) when not part of a workflow. You can also tap **REPORTS** on the initial app screen to view reports for a patient without having to start a session.

Note: Any report generated in Demo mode is labeled **DEMO MODE** at the top of the report.

Using the Reports screen

1. Select a report type from the drop-down list. The default is Session Report.
2. (Optional) Use the filter options to find the patient data you want.
3. Tap the desired patient session record in the scrollable list. The most recent patient session record appears at the top.

4. Tap one of the buttons shown in Table 13.

Table 13. Report buttons.

Button	Description
	View button – to view the report in PDF format. This action takes you out of the clinician programmer app. The app runs in the background and your session data remains intact.
	Share button – to send the PDF report to another destination. This action takes you out of the clinician programmer app. The app runs in the background and your session data remains intact.
	Download button – to download and save the PDF report on the clinician tablet. This action saves a PDF file in the default reports folder on the clinician tablet.
	Delete button – to delete the selected patient session record from the clinician tablet. When a session record is deleted, the data no longer exists on the clinician tablet. However, deleting a session record does not delete any existing PDF files created from that record and stored on the clinician tablet.

Patient Data Service app

The Patient Data Service app manages the storage of patient session data in a locally stored database on the clinician tablet.

The Patient Data Service app can be used to access or delete patient session data for any patients whose Medtronic devices have been programmed using any clinician programmer app on the clinician tablet.

The Patient Data Service app can also be used to adjust (or disable) the PDF auto-delete feature, which automatically deletes PDF files from the default reports folder on the clinician tablet after a designated amount of time.

Finding and opening the app

1. Navigate to the Apps on the clinician tablet.
2. Find the icon for the Patient Data Service app.
3. Tap the icon to open the app.

Accessing or deleting a patient session record

1. In the Patient Data Service app, use the filter options to find the patient data you want.
2. Tap the desired patient session record in the scrollable list. The most recent patient session record appears at the top.
3. Tap one of the buttons shown in Table 14.

Table 14. Patient Data Service buttons.

Button	Description
	Import button – to import the selected session record into the Reports screen in the clinician programmer app. Refer to "Using the Reports screen" on page 55.
	Delete button – to delete the selected patient session record from the clinician tablet. When a session record is deleted, the data no longer exists on the clinician tablet. However, deleting a session record does not delete any existing PDF files created from that record and stored on the clinician tablet.

Deleting multiple patient session records at once

1. In the Patient Data Service app, use the filter options to find the patient data you want.
2. Tap and hold on any record in the results list to activate the multiple record selection feature.
3. Tap the applicable checkboxes to select all records you want to delete from the scrollable list.
4. Tap the **Trash** icon () to delete the records. The number of records to be deleted is shown next to the icon.

Note: To cancel without deleting records, tap the back arrow in the upper left corner.

Setting the PDF auto-delete feature

Use the auto-delete feature to automatically delete PDF files stored on the clinician tablet based on a timeframe you select. The auto-delete feature is enabled by default with a preset amount of time.

This feature only deletes PDF files from the default reports folder on the clinician tablet. It does not affect patient session data stored in the database.

1. In the Patient Data Service app, tap the **Settings** button (⚙️).
2. Select **PDF Auto Delete**.
3. Select the appropriate auto-delete settings, and then tap **UPDATE**.

Viewing System Information

The **About System** screen is accessed from the Side Menu (≡).

This screen provides the following information:

- System component information (such as model numbers, serial numbers, firmware numbers, version numbers) for the device (ENS or INS), clinician tablet, clinician programmer app, controller (patient programmer), recharger, and communicator. Includes the component information for the controller (patient programmer), recharger, and communicator that most recently communicated with the neurostimulator or was paired with the clinician tablet.
- Trademarks and licenses information.

This screen provides the following capabilities:

- Ability to change the device date and time for the ENS and INS.
- Ability to change the INS implant date.
- Ability to generate a Medtronic Data Report.

Viewing information about system components

On the **About System** screen, tap the tab for the component information you want to view.

Viewing information about trademarks and licenses

On the **About System** screen, tap the **Licenses** tab, and then tap the button for the information you want to view.

Changing device date/time or implant date information

The **Device** tab on the **About System** screen shows the device date, device time, and implant date (for INS models). You can change the date and time information if needed.

- Tap the date that you want to change. Use the calendar selector to select a new date, and then tap outside the calendar.
- Tap the time that you want to change. Use the time selector to select a new time, and then tap outside the time.

Note: Changing the implant date does not impact the clock inside rechargeable neurostimulators that determines ERI and EOS.

Compiling a Medtronic Data Report

The **Medtronic Data Report** tab on the **About System** screen is used to compile a report containing information about the Medtronic data collected from the neurostimulator.

- Tap the **COMPILE REPORT** button (D) to compile the Medtronic Data Report. You cannot cancel the action after it is initiated.

To access previously compiled Medtronic Data Reports, refer to "Working with Reports" on page 55.

Viewing Alerts

The clinician programmer app displays alerts to inform you of specific events or changes in the system. There are three levels of alerts in the clinician programmer app: warning, caution, and information.

When alerts are present, an **ALERTS** button displays on the screen. The button shows the number of alerts that need to be acknowledged. Read and follow the instructions given in the alert.

Working with alerts

- To view alerts, tap the **ALERTS** button in the lower left corner of the screen. Each alert appears in a separate pop-up box.
 - Swipe left or right to view individual alerts.
- To close an alert without dismissing it, tap the **X** (or anywhere outside of the alert pop-up box).
- To dismiss an alert, tap the **DISMISS** button. The alert will be marked as a dismissed alert.
- When all alerts are dismissed, the **ALERTS** button reads "NO ALERTS". Those alerts, however, now exist as dismissed alerts and can still be viewed.
- To view dismissed alerts, tap the **NO ALERTS** button, and then tap the **SHOW DISMISSED** button.
- To restore a dismissed alert, tap the **RESTORE** button.

Troubleshooting

Approach troubleshooting conservatively. Prior to performing invasive procedures, ensure that all noninvasive solutions have been considered. If additional assistance is needed, contact Medtronic using the contact information listed on the back cover of this manual.

The clinician programmer app displays messages to inform you of issues in the system. Follow the instructions given in the message.

Data validation error message

During initial interrogation, the app performs a data validation check to ensure that the data in the neurostimulator is valid and can be understood. If data is found to be invalid, the app displays a data validation error message and provides the option to clear the invalid data.

Power on reset (POR) message

If a power on reset message displays at initial interrogation, this indicates that one or more POR events occurred since the last programming session. If a power on reset message displays while in a current programming session, the session will end and any unsaved data will be lost.

Out of regulation (OOR) message

An out of regulation (OOR) message indicates that the delivered amplitude may be lower than the programmed amplitude for the indicated programs because of low impedance. If an out of regulation situation occurs, decrease settings or reduce the number of electrodes for the indicated programs.

Overcharge message

An overcharge message indicates that the battery in a rechargeable neurostimulator has been overcharged. Upon confirmation of the message, the session will continue.

Overdischarge message

An overdischarge message indicates that the patient allowed the battery in a rechargeable neurostimulator to completely discharge, and telemetry is no longer available from the neurostimulator. Refer to "Physician Recharge mode" on page 64 for instructions on restoring the device.

Rapid discharge message

A rapid discharge message indicates that the rechargeable neurostimulator needs to be recharged at least once a day to provide 24-hour stimulation. This issue may be due to high parameter settings or reduced battery capacity. Parameter adjustment or neurostimulator replacement may be needed to increase therapy time between charges.

Elective Replacement Indicator (ERI) message

When an implanted neurostimulator reaches its recommended replacement date, a message displays. Refer to the *System Eligibility, Battery Longevity, Specifications reference manual* for information about neurostimulator battery longevity.

End Of Service (EOS) message

When an implanted neurostimulator reaches its end of service, a message displays. At EOS, therapy is not available and the device

must be replaced. Refer to the *System Eligibility, Battery Longevity, Specifications reference manual* for information about neurostimulator battery longevity.

Using impedance measurements for troubleshooting

Refer to "Checking System Performance" on page 39.

Communicator and clinician tablet troubleshooting

See Table 15 on page 62 for a list of potential issues that can occur with the communicator or the clinician tablet, as well as possible solutions for resolving those issues.

Table 15. Communicator and clinician tablet troubleshooting.

Problem	Possible Solutions
The communicator cannot communicate with the neurostimulator.	■ The communicator may not be positioned correctly over the neurostimulator. – Position the target symbol (⊕) on the back of the communicator so that it is centered over and facing the neurostimulator. ■ The communicator may be too far away from the neurostimulator. – Move the communicator closer to the neurostimulator. ■ Metal surfaces can interfere with communication between the communicator and the neurostimulator. – If the communicator is on a metal table or a metal tray, move the communicator to a nonmetal surface. ■ There could be radio-frequency (RF) interference. – Move the communicator closer to the neurostimulator. ■ The communicator may have interrogated a different device. – Place the communicator directly over the neurostimulator (use a sterile barrier, if applicable) and tap FIND DEVICE on the clinician programmer app.
Cannot pair the communicator with the clinician tablet.	When attempting to connect through BLUETOOTH® wireless technology: ■ Make sure the communicator and clinician tablet are within range of each other (within 3 meters [10 feet]). ■ Make sure that the BLUETOOTH® wireless technology is enabled on the clinician tablet. Refer to "Communicator icons on the tablet status bar" on page 19 for what the communicator icons mean when displayed in the tablet status bar.
The communicator cannot communicate wirelessly with the clinician tablet during a programming session.	Possible reasons: ■ You are in an environment where multiple devices are using BLUETOOTH® wireless technology and thereby creating interference. ■ You have moved into an environment where the BLUETOOTH® wireless technology is prohibited. Possible solutions: ■ Use the USB connector cable to connect the communicator to the clinician tablet. ■ Make sure that the BLUETOOTH® wireless technology is enabled on the clinician tablet. Refer to "Communicator icons on the tablet status bar" on page 19 for what the communicator icons mean when displayed in the tablet status bar.

Problem	Possible Solutions
The communicator battery level is low.	■ Open the battery case and replace the batteries. – If needed, refer to the instructions for replacing the communicator batteries in the <i>Model 8880T2 Communicator technical manual</i>. ■ After replacing the communicator batteries, the clinician programmer app reconnects with the communicator using the BLUETOOTH® wireless technology and resumes the session.
The clinician tablet, the clinician programmer app, or the communicator is unresponsive.	■ Turn off the power for the clinician tablet or the communicator, then turn the power on. ■ If you continue to have communication problems or cannot use the app to program, contact Medtronic using the contact information listed on the back cover of this manual.
The communicator is damaged or overheats.	Use a different communicator: ■ Ensure that the damaged communicator is turned off before using another communicator. ■ If using a replacement communicator prior to starting a programming session, use the USB connector cable and pair the replacement communicator to the clinician tablet. – Refer to "Pairing the communicator to the clinician tablet" on page 12. - Or – If replacing the communicator while in a programming session, use the USB connector cable to connect the communicator to the clinician tablet and continue with the programming session.
If the communicator requires repair, is damaged, or is nonfunctional and a replacement is needed, contact Medtronic using the contact information listed on the back cover of this manual.	

Physician Recharge mode

Physician Recharge mode is only applicable to neurostimulators that are rechargeable (37711, 37712, 37713, 37714, 97712, 97713, 97714).

If the patient allows the implanted neurostimulator to overdischarge, telemetry is no longer available from the neurostimulator. The patient will not be able to charge the neurostimulator battery using the charging system.

The patient must return to the clinic, and the physician must charge the neurostimulator battery using the Physician Recharge mode on the recharger. Typically, the neurostimulator will return to normal charging mode in less than 60 minutes. If it does not, the process should be repeated. If the neurostimulator cannot be charged, it must be replaced.

Notes:

- If the battery is successfully recharged, it regains functionality but may require more frequent charging sessions by the patient.
- The neurostimulator battery can be restored with the Physician Recharge mode only twice. If the battery is allowed to overdischarge a third time, the neurostimulator will reach end of service.

Performing a physician charging session

1. Palpate the patient's skin to determine the best location for the charging antenna, and place the antenna over the neurostimulator.
2. Simultaneously press and hold down the **Start** key and **Audio** key located on the front of the recharger (Figure 5) for 10

seconds. A screen with a 60-minute timer will appear, and charging begins.

Figure 5. Recharger.

3. Continue charging until the recharger displays a new screen indicating the battery charge level and charging efficiency (Figure 6).


Figure 6. Screen displaying charge level and charging efficiency.

4. Repeat steps 1-3 if the new screen does not appear within 60 minutes.
5. Use the clinician programmer to check settings and to turn stimulation on.

Patient programmer messages

The patient programmer displays text and iconic error and informational messages. The applicable messages vary based on neurostimulator model. Some messages provide an error code and tell the patients to contact their physician. These error codes and their troubleshooting procedures are displayed in Table 16.

Table 16. Contact Physician error codes on the patient programmer.

Error code	Explanation
536	Moisture was detected inside the ENS. The ENS must be dried with a towel, the battery removed, and the battery compartment allowed to air dry at room temperature.
537	The Emergency stop key on the ENS is active. Check the key and replace the device if needed.
574	No programs or groups were saved by the clinician programmer. The neurostimulator must be reprogrammed.
00R	Out of regulation situation has occurred. Decrease settings or reduce the number of electrodes.
POR	The neurostimulator has undergone a power on reset. No therapy is available until the device is re-activated.
556/589	The INS battery was overcharged or the wrong batteries were placed in the ENS. The patient must come into the clinic, and the clinician programmer must be used to determine device status.
575/578	Invalid settings were detected in the INS/ENS. The device must be reset and valid settings entered.

Error code	Explanation
0 to 250	Invalid settings were detected in the patient programmer. The patient should remove the programmer batteries and re-insert them after a few seconds. The error code should disappear. If it reappears, the patient programmer may need to be replaced.

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© Medtronic 2022
All Rights Reserved