Medtronic**Cardioblate™ BP2 49321**

Surgical Ablation Device
 Устройство за хирургична аблация
 Pístroj pro chirurgickou ablaci
 Kirurgisk ablationsenhet
 Chirurgische Ablationsklemme
 Συσκευή χειρουργικής κατάλυσης
 Dispositivo de ablación quirúrgica
 Kirurgiline ablaatsiooniseade
 Kirurginen ablaatiolaite
 Appareil d'ablation chirurgicale
 Uređaj za kiruršku ablaciju
 Sebészi ablációs eszköz
 Dispositivo per ablazione chirurgica
 Chirurginės abliacijos prietaisas
 Kirurğiskās ablācijas ierīce
 Помагало за хируршка аблација
 Chirurgisch ablatie-apparaat
 Kirurgisk ablasjonsutstyr
 Urządzenie do ablacji chirurgicznej
 Dispositivo de ablação cirúrgica
 Dispositivo de ablação cirúrgica
 Dispozitiv de ablație chirurgicală
 Устройство для хирургической аблации
 Zariadenie na chirurgickú abláciu
 Pripomoček za kirurško ablacijo
 Uređaj za hiruršku ablaciju
 Enhet för kirurgisk ablation
 Cerrahi Ablasyon Cihazı
 Пристрій для хірургічної абляції

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49321

2 Product description


The Cardioblate BP2 surgical ablation device is a hand-held, single-use, bipolar, radiofrequency (RF) ablation device for use in cardiac surgery. It has a saline irrigation system to deliver fluid at the contact point between the tissue and electrode to cool the tissue during RF energy delivery. The device is intended for intermittent operation.

The Cardioblate BP2 surgical ablation device (Model 49321) is intended for use with the Covidien™ Valleylab™ FT10 FT series energy platform with TissueFect™ technology (Model VLFT10GEN) with software version 3.0.0 or higher.

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Rated accessory voltage: 382 V_p

Figure 1.

- 1 Adjustable parallel jaws, 7 cm (2.76 in) in length
- 2 7 cm (2.76 in) electrodes embedded in a porous polymer, 5 cm (1.97 in) of which is malleable
- 3 300° rotation of jaw assembly
- 4 Malleable neck
- 5 Trigger
- 6 Trigger lock with raised thumb slide
- 7 Molded handle

- 8 Saline fluid delivery tubing 304.8 cm (120 in) in length with a female luer for connection to standard IV tubing systems
 - 9 Electrical cable 304.8 cm (120 in) in length that terminates into a multi-pin connector
-

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2. Open the package and transfer the device onto the sterile field using aseptic technique.

3. Carefully uncoil the saline delivery tubing and the electrical cable to avoid kinking and tangling.


Note: 0.9% normal saline is required because it provides a low impedance pathway for the radiofrequency energy to the tissue.

4. If only a bipolar device is being used, attach the female luer connector of the saline delivery tubing to an IV tubing set and a 1000 mL bag of sterile 0.9% normal saline.

5. If both a bipolar device and a monopolar device are being used, perform the following steps:


- a. Remove the sterile caps from a "Y" connector (provided with the Cardioblate BP2 system, Model 49322) (Figure 2).
- b. Attach the female luer connectors of each saline delivery line to the "Y" connector.
- c. Connect the "Y" connector to an IV tubing set and a 1000 mL bag of sterile 0.9% normal saline.

Figure 2.

6. Insert the 1000 mL bag of 0.9% normal saline into a commercially available pressure cuff and inflate the cuff to between 150 and 300 mmHg.
7. Open the IV line and squeeze the trigger to purge the system of air until fluid is seen coming from the electrodes on the bipolar device. Release the trigger to stop saline flow.

Figure 3.

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Table 1. Compatible generator specifications

Model	VLFT10GEN
RF power output (bipolar mode)	15 to 40 W
RF power output frequency	434 ±10% kHz continuous sinusoid

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8. With the Medtronic logo on the connector facing up, firmly insert the multi-pin connector into the LigaSure™/Bipolar receptacle on the energy platform front panel.
9. Shape the electrodes as required for the targeted tissue. Engage the trigger lock and squeeze the trigger to bring the electrodes together. Grip the electrodes with both hands to shape them at the same time. Ensure that the electrodes are aligned with no more than 4 mm distance between them.
10. Release the trigger lock by moving it backward and open the electrodes. Use the malleable neck and 300° roll to assist in positioning the electrodes on the targeted tissue.
11. Once the electrodes are correctly positioned, engage the trigger lock and squeeze the trigger to bring the electrodes in full contact with the tissue. The lock must be engaged to ensure accurate transmural readings.
12. Confirm that the energy platform touchscreen indicates that the Cardioblate bipolar device is selected for use.
13. Position the device electrical cable to avoid contact with the patient or other leads.
14. Before tissue ablation, squeeze the trigger and ensure that saline is flowing freely from the electrodes.

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Note: Suctioning of excessive fluid at the ablation site is recommended to increase the effectiveness of RF ablation.

Note: In bipolar and enhanced bipolar modes, the energy platform monitors tissue impedance and adjusts the amount of power delivered to the tissue between the active electrodes to create transmural lesions.

Note: Refer to the energy platform user guide for information on high impedance and low impedance alert situations.
15. Activate energy delivery in the device using the **Start RF** button or the foot pedal provided with the energy platform. The energy flows between the electrodes.

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Tissue thickness (7 cm length)	Energy delivered
1 mm	477 J ± 21 J
3 mm	511 J ± 41 J
5 mm	543 J ± 41 J

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16. The achievement of a transmural lesion in the bipolar mode is indicated by the energy platform tone changing from a beeping to a steady tone. The green transmural symbol will appear on the screen next to the power output. RF energy will continue until the surgeon stops ablating or the maximum time-out is reached.

Note: To prevent the porous jaws from drying out and clogging, maintain the jaws immersed in saline or covered with a wet sponge between ablations.
17. When the surgical procedure is complete, disconnect the device from the energy platform and discard the device in accordance with local regulations and hospital procedures.

9 The following disclaimer of warranty applies to customers outside the United States:

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ALTHOUGH THE CARIOBLATE BP2 SURGICAL ABLATION DEVICE, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION

SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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Nicholas
Madson

Digitally signed by Nicholas
Madson
Date: 2021.11.29 12:45:14 -06'00'

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
2017/MAY/29 at 3:30 p.m. Doc number: M972539A001 Rev. 1A [global]
Printspec e - 3.75 x 8.5 inches

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USA Caution: Federal law (USA) restricts
this device to sale by or on the order of a
physician.

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60831

Description

The Cardioblate BP2 surgical ablation device is a hand-held, single-use, bipolar, radiofrequency ablation device for use in cardiac surgery. It has a saline irrigation system to deliver fluid at the contact point between the tissue and electrode to cool the tissue during radiofrequency energy delivery. The device is intended for intermittent operation.

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2. Indications for Use

The Cardioblate BP2 surgical ablation device is intended to create lesions in cardiac tissue during cardiac surgery using radiofrequency energy for the treatment of cardiac arrhythmias. Any supraventricular arrhythmia (such as atrial fibrillation, atrial flutter, or atrial tachycardia) that may be resolved by the application of endocardial or epicardial radiofrequency energy can be treated by the Cardioblate BP2 surgical ablation device.

3. Contraindications

The Cardioblate BP2 surgical ablation device should not be used for:

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5. Warnings

Do not activate device if saline is not flowing freely. The cooling of the outer tissue prevents the charring and tissue desiccation that may occur when applying radiofrequency energy without irrigation.

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Use caution to avoid trauma to tissues not within the target area of ablation.

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Separation of electrodes from tissue contact during activation will result in a high impedance shutdown placing generator power condition to OFF. Placing electrodes in direct contact during radiofrequency energy activation results in a low impedance shut down and places the generator power condition to OFF. In both situations, three short audible alarm beeps will be heard. Refer to the Cardioblate® generator technical manual for reset instructions.

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7. Connections to External Devices

1. Connection to the Cardioblate generator is as follows:
 - a. Position the multi-pin connector of the bipolar device so the alignment indicator matches the generator port indicator.

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Note: Connect to RF generators (only IEC / EN 60601-1 Type CF) with the following specifications:

	Model 60880/60890	Model 68000
RF power output (Transmurality Mode)	25 - 40 W (automatically varying)	15 - 40 W (automatically varying)
RF power output frequency	479 - 489 kHz	468 - 478 kHz
Maximum RF voltage, between electrodes (Rated Accessory Voltage)	159 Vrms	180 Vrms $\pm 10\%$

Note: The rated accessory voltage is limited by the Cardioblate generator.

Note: If additional devices are to be used, refer to their specific instructions for connections to external devices.

2. Sterile normal saline, commercially available pressure cuff, and IV tubing set are required for use with the bipolar device. These are not included with the Cardioblate BP2 surgical ablation device.

8. Instructions for Use

1. Inspect the package and device to ensure the expiration date has not passed and damage has not occurred to the device during shipping and handling. Do not use the device if the packaging or device is damaged or if the expiration date has elapsed. Open the package and transfer the device onto the sterile field utilizing aseptic technique.
2. Uncoil both the saline delivery tubing and the electrical cable carefully to avoid kinking and tangling.
3. Attach the saline delivery tubing female luer connector to an IV tubing set and 1000 mL bag of sterile normal saline. (0.9% normal saline is recommended because it provides a low impedance pathway for the radiofrequency energy to the tissue.)
4. If both a bipolar and monopolar device are being used, remove the sterile caps from the "Y" connector (Figure 2) provided and attach the female luer of both devices. Then insert the IV line into the "Y" connector.


Figure 2.

5. Insert a 1000 mL bag of 0.9% normal saline into commercially available pressure cuff and inflate to between 150 and 300 mm Hg.
6. Open the IV line and squeeze the trigger to purge the system of air until fluid is seen coming from the electrodes on the bipolar device. Release the trigger to stop IV flow.


Figure 3

7.
 - a. Attach the bipolar electrical cable to the generator per instructions outlined in "Connections to External Devices."
 - b. Shape the electrodes as required for the targeted tissue. Engage the lock and squeeze the trigger to bring the electrodes together. Grip the electrodes with both hands to shape them. Care should be taken to ensure the electrodes are as parallel as possible. If the electrodes are out of alignment greater than 4 mm, the ability of the device to create a transmural lesion across the entire electrode length may be affected.
 - c. Release the lock by moving it backward and open the electrodes. Use the flexible neck and 300" roll to assist in positioning the electrodes on the targeted tissue. Once the electrodes are correctly positioned, engage the lock and squeeze the trigger to bring the electrodes in full contact with the tissue. The lock must be engaged to ensure accurate transmural readings.
 - d. Confirm the generator selector is placed on the bipolar device.
 - e. Position the device electrical cable to avoid contact with the patient or other leads.
8. Prior to tissue ablation:
 - a. Squeeze the trigger and make certain that saline is flowing freely out of the electrodes. Continuous flow of saline is necessary during ablation to cool the outer tissue. Cooling of the outer tissue prevents charring and tissue desiccation and allows radiofrequency energy to be delivered deeper into the tissue to more effectively create transmural lesions.
Note: While irrigation increases the effectiveness of radiofrequency ablation, ablating in a pool of fluid can have a detrimental effect. Radiofrequency energy can shunt between the two electrodes and the transmural algorithm may be less effective. Suctioning of excessive fluid at the ablation site is recommended to prevent this from occurring.
 - b. When using a bipolar device, the radiofrequency ablation generator will default to the TRANSMURALITY control mode.
 - c. In the TRANSMURALITY mode, all generator parameters will be automatically adjusted.
9. The current is activated in the device using the RF energy button or the foot switch provided with the Cardioblate generator. The current flows between the electrodes.

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Tissue Thickness (7 cm length)	Energy Delivered (W)
1 mm	439±40
3 mm	505±25
5 mm	557±28

10. The lesion becomes visible by a white coloration on the tissue. In TRANSMURALITY control mode, the generator will emit a solid tone once a transmural lesion is created. This is the indication to the surgeon to stop ablating by releasing the foot switch or pressing the RF energy power button on the generator.


Figure 4. Typical Tissue Impedance Values with Corresponding Power Levels. (Approximate total energy delivered 570 W)

1. Average Impedance (Ohms)
 2. Time (msec)
 3. Power (Watts)
11. Ablating beyond recommended ablation times may result in tissue perforation.
Note: To prevent the porous jaws from drying out and clogging, maintain the jaws immersed in saline or covered with a wet sponge between ablations.
12. Upon completion of the surgical procedure, disconnect the device from the generator and discard.

9. The following disclaimer of warranty applies to United States customers only:

Disclaimer of Warranty

ALTHOUGH THE CARDIABLATE BP2 SURGICAL ABLATION DEVICE, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY

USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

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10. The following disclaimer of warranty applies to customers outside the United States:

Disclaimer of Warranty

ALTHOUGH THE CARDIOBLATE BP2 SURGICAL ABLATION DEVICE, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

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J. Kasper Smith 24 Jun 2021
 JENNIFER KASPER SMITH
 NOTARY PUBLIC - MINNESOTA
 MY COMMISSION EXPIRES 01/31/25

As a registered Notary Public for the State of Minnesota,
 I hereby certify that this is a true and accurate copy of
 the original document. The original is filed at
 Medtronic, Inc.

Instructions for Use • Инструкции за употреба • Návod k použití •
 Brugsanvisning • Gebrauchsanweisung • Οδηγίες χρήσης •
 Instrucciones de uso • Kasutusjuhend • Käyttöohjeet • Mode d'emploi •
 Upute za upotrebu • Használati útmutató • Istruzioni per l'uso •
 Naudojimo instrukcijos • Lietošanas pamācība • Упатство за употреба
 • Gebruiksaanwijzing • Bruksanvisning • Instrukcja użytkowania •
 Instruções de utilização • Instruções de utilização • Instrucțiuni de
 utilizare • Инструкция по эксплуатации • Pokyny na používanie •
 Navodila za uporabo • Uputstva za upotrebu • Bruksanvisning •
 Kullanım Talimatları • Інструкція з використання

CE
 0123

The following list includes trademarks or registered trademarks of Medtronic in the United States and possibly in other countries. All other trademarks are the property of their respective owners.

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49351

2 Product description

The Cardioblade Gemini-s surgical ablation device is a hand-held, single-use, bipolar, radiofrequency (RF) ablation device for use in cardiac surgery. It has a saline irrigation system to deliver fluid at the contact point between the tissue and electrode to cool the tissue during RF energy delivery. The device is intended for intermittent operation.

The Cardioblade Gemini-s surgical ablation device (Model 49351) is intended for use with the Covidien™ Valleylab™ FT10 FT series energy platform with TissueFect™ technology (Model VLFT10GEN) with software version 3.0.0 or higher.

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Rated accessory voltage: 382 V_p

3 Indications for use

The Cardioblade Gemini-s surgical ablation device is intended to create lesions in cardiac tissue during cardiac surgery using radiofrequency energy for the treatment of cardiac arrhythmias. Any supraventricular arrhythmia (such as atrial fibrillation, atrial flutter, or atrial tachycardia) that may be resolved by the application of endocardial or epicardial radiofrequency energy can be treated by the Cardioblade Gemini-s surgical ablation device. The system is indicated for use, under direct or endoscopic visualization, in surgical procedures, including minimally invasive surgical procedures.

4 Contraindications

The Cardioblade Gemini-s surgical ablation device should not be used for:

- Patients that have active endocarditis at time of surgery
- Ablation in a pool of blood (for example, through a purse string suture on a beating heart) (Effects of this type of ablation are unknown.)

5 Adverse effects

The following known adverse effects are associated with the use of the product:

- Abrasion
- Arrhythmia or EKG/ECG changes
- Burn
- Electrical shock
- Infection
- Ischemia
- Organ dysfunction (cardiac)
- Perforation
- Tissue damage

The following are possible adverse effects related to the ablation of cardiac tissue in combination with open heart surgery:

- Atrial lead dislodgement
- Atrioesophageal fistula
- Bleeding
- Cardiac perforation/tamponade
- Conduction disturbances (SA or AV node)
- Esophageal injury
- Extension of extracorporeal bypass
- Major complication (death)
- Mediastinitis
- Myocardial infarction in the context of cardiac ablation
- Pericarditis
- Phrenic nerve paralysis
- Pleural effusion
- Pulmonary vein stenosis
- Stroke or transient ischemic attack (TIA) post-ablation
- Vagal nerve injury
- Valve leaflet damage

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- Damage to the physiological conduction system
- Peripheral tissue ablation due to breach of insulation covering

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Use the device with caution in the presence of pacemakers or other active implanted devices. Interference produced by electrosurgical devices can affect the operation of the active implanted device or cause permanent damage to the active implanted device. Consult the active implanted device manufacturer or the responsible hospital department before using electrosurgical devices in patients with active implanted devices.

7 Precautions

Use caution to visualize the intended area of ablation to avoid trauma to tissues not within the intended area of ablation.

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Applicable to products with phthalates identified on the product label:
This product contains phthalates.

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Warning: Failure to follow instructions may result in tissue damage and insufficient ablation.

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2. Open the package and transfer the device onto the sterile field using aseptic technique.
3. Carefully uncoil the saline delivery tubing and the electrical cable to avoid kinking and tangling.
Note: 0.9% normal saline is required because it provides a low impedance pathway for the RF energy to the tissue.
4. Attach the female luer connector of the saline delivery tubing to an IV tubing set and a 1000 mL bag of sterile 0.9% normal saline.
5. Insert the 1000 mL bag of 0.9% normal saline into a commercially available pressure cuff and inflate the cuff to between 150 and 300 mmHg.
6. Open the IV line and squeeze the handle to purge the system of air until fluid is seen coming from the electrodes on the device. Release the handle to stop saline flow.

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Table 1. Compatible generator specifications

Model	VLFT10GEN
RF power output (enhanced bipolar mode)	15 to 55 W
RF power output frequency	434 ±10% kHz continuous sinusoid

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7. With the Medtronic logo on the connector facing up, firmly insert the multi-pin connector into the **LigaSure™/Bipolar** receptacle on the energy platform front panel.

8. Connect the guides to the device (Figure 3) to assist in positioning the electrodes on the targeted tissue.
9. Confirm that the energy platform touchscreen indicates that the Cardioblate bipolar device is selected for use.
10. Position the device electrical cable to avoid contact with the patient or other leads.
11. Before tissue ablation, squeeze the handle and ensure that saline is flowing freely from the electrodes.
Warning: Do not activate the device if saline is not flowing freely. Applying radiofrequency energy without irrigation can result in charring, tissue desiccation, or insufficient ablation.
Note: Suctioning of excessive fluid at the ablation site is recommended to increase the effectiveness of RF ablation.
Note: In bipolar and enhanced bipolar modes, the energy platform monitors tissue impedance and adjusts the amount of power delivered to the tissue between the active electrodes to create transmural lesions.
Note: Refer to the energy platform user guide for information on high impedance and low impedance alert situations.
12. Once the electrodes are correctly positioned, squeeze the handle until the parallel jaws lock and the electrodes are in full contact with the tissue. The locking mechanism built into the handle must be fully engaged to ensure accurate transmural readings.
13. Activate energy delivery in the device by either pressing the **Start RF** button, squeezing the RF trigger on the device, or using the foot pedal provided with the energy platform. The energy flows between the electrodes. When using the RF trigger, the energy flows between the electrodes only when the handle is in the locked position.

Table 2. Anticipated energy delivery at various tissue thicknesses

Tissue thickness (6 cm length)	Energy delivered
2 mm	696 J $\pm$ 109 J
6 mm	784 J $\pm$ 41 J
10 mm	825 J $\pm$ 74 J

Note: The energy delivery duration is automatically limited by the energy platform.

Note: Tissue thicknesses of 10 mm or greater may require additional ablation cycles to ensure transmural across the entire length of the electrodes.

14. The achievement of a transmural lesion in the bipolar mode is indicated by the energy platform tone changing from a beeping to a steady tone. The green transmural symbol will appear on the screen next to the power output. RF energy will continue until the surgeon stops ablating or the maximum time-out is reached.

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15. After tissue ablation, perform the following steps:
 - a. Press the handle lock release on the handle to unlock the parallel jaws and allow the device to be retracted from the intended tissue.
 - b. Disconnect the guides from the device (Figure 4).
16. When the surgical procedure is complete, disconnect the device from the energy platform and discard the device in accordance with local regulations and hospital procedures.

Table 3 and Table 4 provide examples of lesion depth, width, and length using the Cardioblate Gemini-s surgical ablation device:

Table 3. In vitro lesion evaluation study

	Cardiac Tissue (n=42)
Max transmural depth (mm)	5.6 $\pm$ 1.2
Lesion length (cm)	4.1 $\pm$ 0.6
Lesion width (mm)	4.5 $\pm$ 0.4

Table 4. In vivo lesion observations in a single porcine animal study

	Cardiac Tissue (n=3)
Max transmural depth (mm)	3.0 ±2.5 ^a
Lesion length (cm)	3.8 ±0.6
Lesion width (mm)	4.3 ±1.2

^a Data was obtained in a nonstatistical sample size, in vivo porcine study with 3 lesion burns.

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Nicholas Madson

Nicholas Madson
Associate Regulatory Affairs Specialist

Nicholas
Madson

Digitally signed by
Nicholas Madson
Date: 2021.11.29
12:48:41 -06'00'

2017/APR/03 at 2:21 p.m. Doc number: M972542A001 Rev. 1A [global]
Printspec ø - 3.75 x 8.5 inches


Medtronic

Cardioblate™ Gemini™

Surgical Ablation Device • Устройство за хирургична аблация • Přístroj pro chirurgickou ablací • Kirurgisk ablationsenhed • Chirurgische Ablationsklemme • Συσκευή χειρουργικής κατάλυσης • Dispositivo de ablación quirúrgica • Kirurgiline ablatiooniseade • Kirurginen ablaatiolaitte • Appareil d'ablation chirurgicale • Uređaj za kiruršku ablaciju • Sebészi ablációs eszköz • Dispositivo per ablazione chirurgica • Chirurginės abliacijos prietaisas • Kirurgiskās ablācijas ierīce • Chirurgisch ablatie-apparaat • Kirurgisk ablasjonsutstyr • Urządzenie do ablacji chirurgicznej • Dispositivo de ablação cirúrgica • Dispozitiv de ablație chirurgicală • Устройство для хирургической аблации • Zariadenie na chirurgickú abláciu • Pripomoček za kirurško ablacijo • Uređaj za hiruršku ablaciju • Enhet för kirurgisk ablation • Cerrahi Ablasyon Cihazı

49260

Varina Phil 30 NOV 2021

As a registered Notary Public for the State of Minnesota,
I hereby certify that this is a true and accurate copy of
the original document. The original is filed at
Medtronic, Inc.

Instructions for Use • Инструкции за употреба •
Návod k použití • Brugsanvisning • Gebrauchsanweisung •
Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend •
Käyttöohjeet • Mode d'emploi • Upute za upotrebu •
Használati útmutató • Istruzioni per l'uso •
Naudojimo instrukcijos • Lietošanas pamācība •
Gebruiksaanwijzing • Bruksanvisning •
Instrukcja użytkowania • Instruções de utilização •
Instrucțiuni de utilizare • инструкция по применению •
Pokyny na používanie • Navodila za uporabo •
Uputstva za upotrebu • Bruksanvisning • Kullanım Talimatları

 **Caution:** Federal law (USA) restricts
this device to sale by or on the order of a
physician.

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49260

Description

The Medtronic Cardioblate Gemini surgical ablation device is a hand-held, single-use, bipolar, radiofrequency ablation device for use in cardiac surgery. It has a saline irrigation system to deliver fluid at the contact point between the tissue and electrode to cool the tissue during radiofrequency energy delivery. The device is intended for intermittent operation.

Sterile, nonpyrogenic, disposable, single use only.

2. Indications for Use

The Cardioblate Gemini surgical ablation device is intended to create lesions in cardiac tissue during cardiac surgery using radiofrequency energy for the treatment of cardiac arrhythmias. Any supraventricular arrhythmia (such as atrial fibrillation, atrial flutter, or atrial tachycardia) that may be resolved by the application of endocardial or epicardial radiofrequency energy can be treated by the Cardioblate Gemini surgical ablation device. The system is indicated for use, under direct or endoscopic visualization, in surgical procedures, including minimally invasive surgical procedures.

3. Contraindications

The Cardioblate Gemini surgical ablation device should not be used for:

- patients that have active endocarditis at time of surgery
- ablation in a pool of blood (eg, through a purse string suture on a beating heart) (Effects of this type of ablation are unknown.)

Adverse Effects

Possible adverse effects related to the ablation of cardiac tissue in combination with open heart surgery are:

- tissue perforation
- extension of extracorporeal bypass
- perioperative heart rhythm disturbances (atrial and/or ventricular)
- postoperative embolic complications
- pericardial effusion or tamponade
- injury to the great vessels
- valve leaflet damage
- conduction disturbances (SA/AV node)
- acute ischemic myocardial event
- thrombus formation
- nerve damage
- unintentional burns
- pericarditis
- pleural effusion
- esophageal perforation
- death
- coronary sinus perforation
- coronary artery spasm
- atrial lead dislodgement
- hypotension
- pulmonary vein stenosis
- cerebrovascular accident
- transient ischemic attack
- blood loss

4. Warnings

Do not activate device if saline is not flowing freely. Cooling of the outer tissue prevents the charring and tissue desiccation that may occur when applying radiofrequency energy without irrigation.

The device may be damaged before or during the procedure by improper handling or other intervening acts. If damaged, the device may fail to function properly, and may result in the following medical complications including, but not limited to:

- damage to the physiological conduction system
- peripheral tissue ablation due to breach of insulation covering

Device shall be removed from the patient before defibrillation.

No representation or warranty is made that failure or cessation of function of the device will not result in an adverse event, or that medical complications (including perforation of cardiac tissue) will not follow the procedure, or that the use of the device will in all cases restore adequate cardiac function.

5. Precautions

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Proper surgical procedures and techniques are necessarily the responsibility of the medical professional. Each surgeon must evaluate the appropriateness of any procedure based on his or her own medical training and experience, and the type of surgical procedure.

Patient and procedure selection is solely a medical responsibility and the outcome is dependent on many variables including patient pathology, surgical and perfusion procedures.

Use caution to visualize the target area to avoid trauma to tissues not within the target area of ablation.

Purge air from system to ensure saline fluid flows from the electrodes during use.

Separation of electrodes from tissue contact during activation will result in a high impedance shutdown placing generator power condition to OFF. Placing electrodes in direct contact during radiofrequency energy activation results in a low impedance shutdown and places the generator power condition to OFF. In both situations, three short audible alarm beeps will be heard. Refer to the Cardioblate 68000 generator technical manual for reset instructions.

Maximum storage temperature: 25°C (77°F)

Storage relative humidity: 5% to 85%, noncondensing

This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death. Multiple uses can result in the occlusion of the irrigation openings and affect the device's performance. To prevent this from occurring, the generator has the capability to detect previously used ablation devices and will not allow their reuse.

It is the responsibility of the user to dispose of the devices in accordance with local regulations and hospital procedures.

Due to the presence of phthalates in the product, the clinician must weigh the medical benefits of product use against the drawbacks of phthalate exposure for male children and pregnant or nursing women.

6. Connections to External Devices

1. Connection to the Cardioblate 68000 Generator is as follows:
 - a. Position the multi-pin connector of the bipolar device so the alignment indicator matches the generator port indicator.

- b. Insert the multi-pin connector of the bipolar device into the device port on the generator until it locks in place. Do not force the connection.

Note: Connect to RF generators (only IEC/EN 60601-1 Type CF) with the following specifications:

	Model 68000
RF power output (Transmurality Mode)	<55 W
RF power output frequency	468 - 478 kHz
Maximum RF voltage, between electrodes (Rated Accessory Voltage)	180 Vrms $\pm$ 10%

Note: If additional devices are to be used, refer to their specific Instructions for Use.

2. Sterile normal saline, commercially available pressure cuff, and IV tubing set are required for use with the bipolar device. These are not included with the Cardioblate Gemini surgical ablation device.

7. Instructions for Use

1. Inspect the package and device to ensure the expiration date has not passed and damage has not occurred to the device during shipping and handling. Do not use the device if the packaging or device is damaged or if the expiration date has elapsed. Open the package and transfer the contents onto the sterile field utilizing aseptic technique.
2. Uncoil both the saline delivery tubing and the electrical cable carefully to avoid kinking and tangling.
3. Attach the saline delivery tubing female luer connector to an IV tubing set and 1000 mL bag of sterile normal saline. (0.9% normal saline is recommended because it provides a low impedance pathway for the radiofrequency energy to the tissue.)
4. Insert a 1000 mL bag of 0.9% normal saline into commercially available pressure cuff and inflate to between 150 and 300 mm Hg.
5. Open the IV line and squeeze the handle to purge the system of air until fluid is seen coming from the electrodes on the Gemini device. Release the handle to stop IV flow.


Figure 2. Jaw/Electrode

6. Attach the Gemini electrical cable to the Cardioblate 68000 Generator per instructions outlined in "Connections to External Devices."
7. Prior to tissue ablation:
 - a. Squeeze the handle and make certain that saline is flowing freely out of the electrodes. Continuous flow of saline is necessary during ablation to cool the outer tissue. Cooling of the outer tissue prevents charring and tissue desiccation and allows radiofrequency energy to be delivered deeper into the tissue to more effectively create transmural lesions.

Note: While irrigation increases the effectiveness of radiofrequency ablation, ablating in a pool of fluid can have a detrimental effect. Radiofrequency energy can shunt between

- the 2 electrodes and the transmural algorithm may be less effective. Suctioning of excessive fluid at the ablation site is recommended to prevent this from occurring.
- b. Confirm the Cardioblate 68000 Generator is set to the bipolar mode.
 - c. Position the device electrical cable to avoid contact with the patient or other leads.
 - d. The mode of operation will default to the TRANSMURALITY control mode.
 - e. In the TRANSMURALITY mode, all generator parameters will be automatically adjusted.
8. Tissue Ablation:
- a. Once the electrodes are correctly positioned, squeeze the handle until the parallel jaws lock and the electrodes are in full contact with the tissue. The locking mechanism built into the handle must be fully engaged to ensure accurate transmural readings.
 - b. The RF energy is activated in the device by either pressing the RF energy button on the Cardioblate 68000 Generator display panel, squeezing the RF Trigger on the device, or using the foot switch provided with the generator. When using the RF Trigger, the current flows between the electrodes only when the handle is in the locked position.

Table 1. Anticipated Energy Delivery at Various Tissue Thicknesses

Tissue Thickness (6 cm length)	Energy Delivered (W)
2 mm	833 ±32
6 mm	889 ±135
10 mm	1023 ±163

- c. The lesion becomes visible by a white coloration on the tissue. In TRANSMURALITY control mode, the Cardioblate 68000 Generator will emit a solid tone once a transmural lesion is created. This is the indication to the surgeon to stop ablating by releasing the foot switch, releasing the RF trigger, or pressing the RF energy button on the generator.

Note: Ablation times will vary due to tissue thickness but will not exceed 45 seconds per ablation cycle.

Note: Ablating beyond recommended ablation times may result in tissue perforation.

Note: To prevent the porous jaws from drying out and clogging, maintain the jaws immersed in saline or covered with a wet sponge between ablations.


Figure 3. Typical Tissue Impedance Values with Corresponding Power Levels. (Approximate total energy delivered 570 W)

1. Average Impedance (Ohms)
 2. Time (msec)
 3. Power (Watts)
 9. After tissue ablation:
 - Press the handle lock release on the handle to unlock the parallel jaws and allow the device to be retracted from the target tissue.
 10. Upon completion of the surgical procedure, disconnect the device from the generator and discard.
- The tables below provide an example of lesion depth, width, and length using the Gemini device:

	In vitro lesion evaluation study	
	Cardiac Tissue (n=42)	Soft Tissue (n=44)
Max Transmurality Depth (mm)	5.6 ±1.2	3.0*
Lesion Length (cm)	4.1 ±0.6	6.0*
Lesion Width (mm)	4.5 ±0.4	2.9 ±1.1

*In vitro transmurality and length variables were controlled by tissue cut to 3 mm thickness and full 6 cm bites, therefore there are no standard deviations on soft tissue lesions.

	In vivo lesion observations in a single porcine animal study			
	Cardiac Tissue (n=3)	Pulmonary Tissue (n=3)	Hepatic Tissue (n=3)	Renal Tissue (n=3)
Max Transmurality Depth (mm)	3.0 ±2.5**	6.5 ±1.9	7.4 ±2.9	8.5 ±2.4
Lesion Length (cm)	3.8 ±0.6	3.8 ±1.0	5.3 ±0.4	4.8 ±1.0
Lesion Width (mm)	4.3 ±1.2	5.7 ±2.0	5.4 ±1.0	3.6 ±1.3

**Data was obtained in a nonstatistical sample size, in vivo porcine study with 3 lesion burns in each tissue. The Gemini device can

achieve transmural of up to 10 mm in cardiac tissue as demonstrated in vitro.

8. The following disclaimer of warranty applies to United States customers only:

Disclaimer of Warranty

ALTHOUGH THE CARIOBLATE GEMINI SURGICAL ABLATION DEVICE, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law, by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

9. The following disclaimer of warranty applies to customers outside the United States:

Disclaimer of Warranty

ALTHOUGH THE CARIOBLATE GEMINI SURGICAL ABLATION DEVICE, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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49313

2 Product description

The Cardioblade surgical ablation pen is a hand-held, single-use, monopolar, radiofrequency (RF) ablation device for use in cardiac surgery. It has a saline irrigation system to deliver fluid at the contact point between the tissue and electrode tip to cool the tissue during RF energy delivery. The device is intended for intermittent operation.

The Cardioblade surgical ablation pen (Model 49313) is intended for use with the Covidien™ Valleylab™ FT10 FT series energy platform with TissueFect™ technology (Model VLFT10GEN) with software version 3.0.0 or higher.

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Rated accessory voltage: 382 V_p

3 Indications for use

The Cardioblade surgical ablation pen is intended to create spot or linear lesions in cardiac tissue during cardiac surgery using RF energy for the treatment of cardiac arrhythmias. Any supraventricular arrhythmia (such as atrial fibrillation, atrial flutter, or atrial tachycardia) that may be resolved by the application of endocardial or epicardial radiofrequency energy can be treated by the Cardioblade surgical ablation pen.

4 Contraindications

The Cardioblade surgical ablation pen should not be used for patients that have active endocarditis at time of surgery.

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The following are possible adverse effects related to the creation of spot or linear lesions in cardiac tissue in combination with open heart surgery: atrial lead dislodgment, atrioesophageal fistula, bleeding, cardiac perforation/tamponade, conduction disturbances (SA or AV node), esophageal injury, extension of extracorporeal bypass, major complication (death), mediastinitis, myocardial infarction in the context of cardiac ablation, pericarditis, phrenic nerve paralysis, pleural effusion, pulmonary vein stenosis, stroke or transient ischemic attack (TIA) post-ablation, vagal nerve injury, and valve leaflet damage.

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4. Connect the patient return electrode to the energy platform according to the instructions in the energy platform user guide and the patient return electrode instructions for use.
5. Carefully uncoil the saline delivery tubing and the electrical cable to avoid kinking and tangling.
Note: 0.9% normal saline is required because it provides a low impedance pathway for the RF energy to the tissue.
6. If only a monopolar device is being used, attach the female luer connector of the saline delivery tubing to an IV tubing set and a 1000 mL bag of sterile 0.9% normal saline.
7. If both a bipolar device and a monopolar device are being used, perform the following steps:
 - a. Remove the sterile caps from a "Y" connector (provided with the Cardioblate BP2 and LP systems, Models 49322 and 49342) (Figure 2).
 - b. Attach the female luer connectors of each saline delivery line to the "Y" connector.
 - c. Connect the "Y" connector to an IV tubing set and a 1000 mL bag of sterile 0.9% normal saline.
8. Insert the 1000 mL bag of 0.9% normal saline into a commercially available pressure cuff and inflate the cuff to between 150 and 300 mmHg.
9. Purge the system of air until fluid is seen coming from the openings on the distal tip of the device (Figure 3). Clamp off the tubing set until ready for use.
Note: Only connect the device to RF generators with the following specifications:

Table 1. Compatible generator specifications

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- Note:** The rated accessory voltage is limited by the energy platform.
10. With the Medtronic logo on the connector facing up, firmly insert the multi-pin connector into the Monopolar 2 receptacle on the energy platform front panel.
 11. Confirm that the energy platform touchscreen indicates that the Cardioblate monopolar device is selected for use.
 12. Position the device electrical cable to avoid contact with the patient or other leads.
 13. Before tissue ablation, open the IV line and ensure that saline is flowing freely from the openings in the distal tip of the device.
Warning: Do not activate the device if saline is not flowing freely. Applying radiofrequency energy without irrigation can result in charring, tissue desiccation, or insufficient ablation.
Warning: Suction excessive fluid at the ablation site before activating the device. Failure to suction excessive fluid can result in unintended burns to the patient.
Caution: The intensity should be set as low as necessary to achieve the desired effect.
Note: Refer to the energy platform user guide for information on high impedance and low impedance alert situations.
Note: Ideal power settings should be determined by the surgeon during tissue ablation.
 14. Set the power setting on the energy platform to between 20 and 30 W.

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Nicholas
Madson

Digitally signed by
Nicholas Madson
Date: 2021.11.29
12:50:28 -06'00'


2017/APR/26 at 8:44 a.m. Doc number: M972543A001 Rev. 1A [global]
Printspec e - 3.75 x 8.5 inches

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
60813

Description

The Cardioblate surgical ablation pen is a hand-held, single-use, monopolar, radiofrequency ablation device for use in cardiac surgery. It has a saline irrigation system to deliver fluid at the contact point between the tissue and electrode tip to cool the tissue during radiofrequency energy delivery. The device is intended for intermittent operation.

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2. Indications for Use

The Cardioblate surgical ablation pen is intended to create spot or linear lesions in cardiac tissue during cardiac surgery using RF energy for the treatment of cardiac arrhythmias. Any supraventricular arrhythmia (such as atrial fibrillation, atrial flutter, or atrial tachycardia) that may be resolved by the application of endocardial or epicardial radiofrequency energy can be treated by the Cardioblate surgical ablation pen.

3. Contraindications

The Cardioblate surgical ablation pen should not be used for patients that have active endocarditis at time of surgery.

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Follow recommended power settings (step 8) for optimal product performance.

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8. Instructions for Use

1. Inspect the package and device to ensure the expiration date has not passed and no damage occurred to the device during shipping and handling. Do not use the device if the packaging or device is damaged or if the expiration date has elapsed. Open the package and transfer the device onto the sterile field utilizing an aseptic technique.
2. Place a large indifferent electrode pad with electrical cable securely on the back of the patient's body. (Refer to generator manual and indifferent electrode instructions for use for specific instructions on electrode placement and skin care.)
3. Attach the electrical cable from the indifferent electrode securely to the appropriate patient grounding connection site on the Cardioblate generator. (Refer to the Cardioblate generator technical manual.)
4. Uncoil both the saline delivery tubing and the electrical cable carefully to avoid kinking and tangling.
5. Attach the saline delivery tubing female luer connector to an IV tubing set and sterile normal saline bag. (0.9% normal saline is recommended because it provides a low impedance pathway for the radiofrequency energy to the tissue.)
If both a bipolar and monopolar device are being used, attach the female luer of both devices to the "Y" connector (Model 10004D) provided. Then insert the IV line into the "Y" connector.
6. Place the sterile normal saline bag into the pressure bag. Inflate the pressure bag to between 150 and 300 mm Hg. Purge the system of air until fluid is seen coming from the openings on the distal tip of the device (Figure 2). Clamp off tubing set until ready for use.


Figure 2.

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Disclaimer of Warranty

ALTHOUGH THE CARDIOLATE SURGICAL ABLATION PEN, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

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Disclaimer of Warranty

ALTHOUGH THE CARDIOBLATE SURGICAL ABLATION PEN, HEREFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

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49314

2 Product description

The Cardioblade XL surgical ablation pen is a hand-held, single-use, monopolar, radiofrequency (RF) ablation device for use in cardiac surgery. It has a saline irrigation system to deliver fluid at the contact point between the tissue and electrode tip to cool the tissue during RF energy delivery. The device is intended for intermittent operation.

The Cardioblade XL surgical ablation pen (Model 49314) is intended for use with the Covidien™ Valleylab™ FT10 FT series energy platform with TissueFect™ technology (Model VLFT10GEN) with software version 3.0.0 or higher.

Sterile, nonpyrogenic, disposable, single use only. Sterilized using ethylene oxide.

Rated accessory voltage: 382 V_p

3 Indications for use

The Cardioblade XL surgical ablation pen is intended to create spot or linear lesions in cardiac tissue during cardiac surgery using RF energy for the treatment of cardiac arrhythmias. Any supraventricular arrhythmia (such as atrial fibrillation, atrial flutter, or atrial tachycardia) that may be resolved by the application of endocardial or epicardial radiofrequency energy can be treated by the Cardioblade XL surgical ablation pen.

4 Contraindications

The Cardioblade XL surgical ablation pen should not be used for patients that have active endocarditis at time of surgery.

5 Adverse effects

The following known adverse effects are associated with the use of the product: abrasion, arrhythmia or EKG/ECG changes, burn, electrical shock, infection, ischemia, organ dysfunction (cardiac), perforation, and tissue damage.

The following are possible adverse effects related to the creation of spot or linear lesions in cardiac tissue in combination with open heart surgery: atrial lead dislodgment, atriopharyngeal fistula, bleeding, cardiac perforation/tamponade, conduction disturbances (SA or AV node), esophageal injury, extension of extracorporeal bypass, major complication (death), mediastinitis, myocardial infarction in the context of cardiac ablation, pericarditis, phrenic nerve paralysis, pleural effusion, pulmonary vein stenosis, stroke or transient ischemic attack (TIA) post-ablation, vagal nerve injury, and valve leaflet damage.

6 Warnings

Do not activate the device if saline is not flowing freely. Applying radiofrequency energy without irrigation can result in charring, tissue desiccation, or insufficient ablation.

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- Tissue "popping" and perforation
- Damage to the physiological conduction system
- Peripheral tissue ablation due to breach of insulation covering

Remove the device from the patient before defibrillation.

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This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death. Multiple uses can result in the occlusion of the irrigation openings and affect the performance of the device. The energy platform has functionality to prevent the reuse of the device.

It is the responsibility of the hospital to ensure that the individuals assigned to operate this device are well-trained in the operation of the device and alert to potential problems.

Do not modify the device.

Do not use the device in the presence of flammable anesthetics, oxidizing gases, or flammable materials used for cleaning or disinfecting, as explosion may occur.

Use the device with caution in the presence of pacemakers or other active implanted devices. Interference produced by electrosurgical devices can affect the operation of the active implanted device or cause permanent damage to the active implanted device. Consult the active implanted device

manufacturer or the responsible hospital department before using electrosurgical devices in patients with active implanted devices.

7 Precautions

Use caution to avoid trauma to tissues not within the intended area of ablation.

Follow recommended power settings for optimal ablation product performance.

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Applicable to products with phthalates identified on the product label:

This product contains phthalates.

Inspect the device and device cable for breaks, cracks, nicks, and other damage before using the device. Do not use the device if any damage is found.

The use of smoke evacuation equipment is recommended. Smoke generated by electrosurgical devices is potentially harmful to patients and users.

8 Instructions for use

Warning: The device is sterile and nonpyrogenic. Inspect each package and device before use. Do not use the device if the use-by date has passed or if the device or package is damaged. Using an expired or damaged device increases the risk of infection.

Warning: Failure to follow instructions may result in tissue damage, charring, tissue desiccation, and insufficient ablation.

Note: Read the energy platform user guide before using the device.

Note: If additional devices are to be used, refer to their specific instructions for use.

Note: Proper surgical procedures and techniques are necessarily the responsibility of the medical professional. Each surgeon must evaluate the appropriateness of any procedure based on his or her own medical training and experience and the type of surgical procedure.

Note: Patient and procedure selection is solely a medical responsibility and the outcome is dependent on many variables including patient pathology, surgical procedures, and perfusion procedures.

Note: Sterile 0.9% normal saline, a commercially available pressure cuff, and an IV tubing set are required to use the device. These items are not included with the device.

1. Inspect the package and the device to ensure that the use-by date has not passed and damage has not occurred during shipping and handling.

Warning: Failure to use aseptic technique increases the risk of contaminating the disposables and harming the patient.

2. Open the package and transfer the device onto the sterile field using aseptic technique.
3. Apply a large patient return electrode (neutral electrode) to the patient according to the instructions in the energy platform user guide and the patient return electrode instructions for use. Choose an application site for the patient return electrode that is as close to the heart as possible and that allows the entire area of the patient return electrode to make secure contact with the patient's skin. The recommended application site is the patient's back.

Note: Only use the following patient return electrodes with the device:

Model	Description
E7507	Valleylab™ REM™ Polyhesive™ Adult Patient Return Electrode
E7507DB	Valleylab REM Polyhesive Adult Patient Return Electrode

Model	Description
E7508	Valleylab REM Polyhesive Adult Cordless Patient Return Electrode
E7509	Valleylab REM Polyhesive Adult Cordless Patient Return Electrode
E7509B	Valleylab REM Polyhesive Adult Cordless Patient Return Electrode

4. Connect the patient return electrode to the energy platform according to the instructions in the energy platform user guide and the patient return electrode instructions for use.
5. Carefully uncoil the saline delivery tubing and the electrical cable to avoid kinking and tangling.

Note: 0.9% normal saline is required because it provides a low impedance pathway for the RF energy to the tissue.

6. If only a monopolar device is being used, attach the female luer connector of the saline delivery tubing to an IV tubing set and a 1000 mL bag of sterile 0.9% normal saline.
7. If both a bipolar device and a monopolar device are being used, perform the following steps:
 - a. Remove the sterile caps from a "Y" connector (provided with the Cardioblate BP2 and LP systems, Models 49322 and 49342) (Figure 2).
 - b. Attach the female luer connectors of each saline delivery line to the "Y" connector.
 - c. Connect the "Y" connector to an IV tubing set and a 1000 mL bag of sterile 0.9% normal saline.
8. Insert the 1000 mL bag of 0.9% normal saline into a commercially available pressure cuff and inflate the cuff to between 150 and 300 mmHg.
9. Purge the system of air until fluid is seen coming from the openings on the distal tip of the device (Figure 3). Clamp off the tubing set until ready for use.

Note: Only connect the device to RF generators with the following specifications:

Table 1. Compatible generator specifications

Model	VLFT10GEN
RF power output (monopolar mode)	1 to 50 W (manually variable)
RF power output frequency	434 ±10% kHz continuous sinusoid

Note: The rated accessory voltage is limited by the energy platform.

10. With the Medtronic logo on the connector facing up, firmly insert the multi-pin connector into the Monopolar 2 receptacle on the energy platform front panel.
11. Confirm that the energy platform touchscreen indicates that the Cardioblate monopolar device is selected for use.
12. Position the device electrical cable to avoid contact with the patient or other leads.
13. Before tissue ablation, open the IV line and ensure that saline is flowing freely from the openings in the distal tip of the device.

Warning: Do not activate the device if saline is not flowing freely. Applying radiofrequency energy without irrigation can result in charring, tissue desiccation, or insufficient ablation.

Warning: Suction excessive fluid at the ablation site before activating the device. Failure to suction excessive fluid can result in unintended burns to the patient.

Caution: The intensity should be set as low as necessary to achieve the desired effect.

Note: Refer to the energy platform user guide for information on high impedance and low impedance alert situations.

Note: Ideal power settings should be determined by the surgeon during tissue ablation.

14. Set the power setting on the energy platform to between 20 and 30 W.

15. Activate energy delivery in the device using the **Start RF** button or the foot pedal provided with the energy platform. Gently drag the distal tip of the device across the intended intracardiac tissue under direct visualization. The energy flows between the tip and the patient return electrode on the patient's body.

Warning: Do not dwell in an intended tissue area too long. This can lead to sudden tissue "popping" and perforation.

16. The lesion becomes visible as a white coloration on the tissue. Move the tip in an oscillating motion at a recommended speed of 1 cm/s. The surgeon determines the depth, width, and location of the ablation.

Note: In vitro testing was conducted using a Valleylab FT10 energy platform (Model VLFT10GEN) and Cardioblate surgical ablation pens to create lesions on bovine myocardium. The testing was performed at 4 combinations of power and time, which resulted in the lesion depth intervals calculated at 95% confidence in Table 2. Results may vary in clinical use due to user technique and cardiac tissue type.

Table 2. Anticipated lesion depth at various power levels

95% confidence interval for 10 s oscillation (1 cm)	95% confidence interval for 20 s oscillation (1 cm)	Power at 300 mmHg pressure on saline irrigation
2.1 to 2.7 mm	3.3 to 4.0 mm	20 W
2.1 to 2.8 mm	5.2 to 6.6 mm	30 W

17. When the surgical procedure is complete, disconnect the device from the energy platform and discard the device in accordance with local regulations and hospital procedures.

9 Disclaimer of warranty

Important: This disclaimer of warranty does not apply in Australia.

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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[Signature]

Nicholas Madson
Nicholas Madson
Associate Regulatory Affairs Specialist

Nicholas
Madson

Digitally signed by
Nicholas Madson
Date: 2021.11.29 12:52:40
+06'00'


2017/JUN/06 at 2:19 p.m. Doc number: M972544A001 Rev. 1A [global]
Printspec e - 3.75 x 8.5 inches

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60814

Description

The Cardioblate XL surgical ablation pen is a hand-held, single-use, monopolar, radiofrequency ablation device for use in cardiac surgery. It has a saline irrigation system to deliver fluid at the contact point between the tissue and electrode tip to cool the tissue during radiofrequency energy delivery. The device is intended for intermittent operation.

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2. Indications for Use

The Cardioblate XL surgical ablation pen is intended to create spot or linear lesions in cardiac tissue during cardiac surgery using RF energy for the treatment of cardiac arrhythmias. Any supraventricular arrhythmia (such as atrial fibrillation, atrial flutter, or atrial tachycardia) that may be resolved by the application of endocardial or epicardial radiofrequency energy can be treated by the Cardioblate XL surgical ablation pen.

3. Contraindications

The Cardioblate XL surgical ablation pen should not be used for patients that have active endocarditis at time of surgery.

4. Adverse Effects

Possible adverse effects related to the creation of spot or linear lesions in cardiac tissue in combination with open heart surgery are:

- tissue perforation
- extension of extracorporeal bypass
- perioperative heart rhythm disturbances (atrial and/or ventricular)
- postoperative embolic complications
- pericardial effusion or tamponade
- injury to the great vessels
- valve leaflet damage
- conduction disturbances (SA/AV node)
- acute ischemic myocardial event
- thrombus formation
- nerve damage
- unintentional burns
- pericarditis
- pleural effusion
- esophageal perforation
- death
- coronary sinus perforation
- coronary artery spasm
- atrial lead dislodgement
- hypotension
- cerebrovascular accident
- transient ischemic attack
- blood loss

5. Warnings

Do not activate ablation if saline is not flowing freely. Cooling of the outer tissue prevents the charring and tissue desiccation that may occur when applying radiofrequency energy without irrigation.

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- tissue "popping" and perforation
- damage to the physiological conduction system
- peripheral tissue ablation due to breach of insulation covering

The device shall be removed from the patient before defibrillation.

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Proper surgical procedures and techniques are the responsibility of the medical professional. Each surgeon must evaluate the appropriateness of any procedure based on their own medical training and experience, and the type of surgical procedure.

Patient and procedure selection is solely a medical responsibility and the outcome is dependent on many variables including patient pathology, and surgical and perfusion procedures.

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Make certain that saline fluid flows from the distal tip ports during use.

Follow recommended power settings (step 8) for optimal ablation product performance.

Separation of tip from tissue contact during activation of RF ablation will result in high impedance shutdown switching the generator power condition to STANDBY. Three short audible alarm beeps from the generator indicate a high impedance occurrence. Refer to the Cardioblate generator technical manual for reset instructions.

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This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death. Multiple uses can result in the occlusion of the irrigation openings and affect the device's performance.

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7. Connections to External Devices

Note: Do not attempt to reuse the ablation device with the Cardioblate generator. The generator has the capability to detect previously used ablation devices and will not allow their reuse.

1. Connection to the Cardioblate generator is as follows:
 - a. Position the multi-pin connector of the pen so the alignment indicator matches the monopolar generator port indicator.
 - b. Insert the multi-pin connector of the pen into the device port on the generator until it locks in place. Do not force the connection.

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	Model 60880/60890	Model 68000
RF power output setpoint/limit	25 - 50 W (manually variable)	25 - 50 W (manually variable)
RF power output frequency	479 - 489 kHz	468 - 478 kHz
Maximum RF voltage, between electrodes (Rated Accessory Voltage)	159 Vrms	180 Vrms $\pm 10\%$

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Note: If additional devices are to be used, refer to their specific directions for connection information.

2. Sterile normal saline, commercially available pressure bag, and IV tubing set are required.

8. Instructions for Use

1. Inspect the package and device to ensure the expiration date has not passed and no damage occurred to the device during shipping and handling. Do not use the device if the packaging or device is damaged or if the expiration date has elapsed. Open the package and transfer the device onto the sterile field utilizing an aseptic technique.
2. Place a large indifferent electrode pad with electrical cable securely on the back of the patient's body (refer to generator manual and indifferent electrode directions for use for specific instructions on electrode placement and skin care).
3. Attach the electrical cable from the indifferent electrode securely to the appropriate patient grounding connection site on the Cardioblate generator (refer to the Cardioblate generator technical manual).
4. Uncoil both the saline delivery tubing and the electrical cable carefully to avoid kinking and tangling.
5. Attach the saline delivery tubing female luer connector to an IV tubing set and sterile normal saline bag. (0.9% normal saline is recommended because it provides a low impedance pathway for the radiofrequency energy to the tissue.)
If both a bipolar and monopolar device are being used, attach the female luer of both devices to the "Y" connector (Model 10004D) provided. Then insert the IV line into the "Y" connector.
6. Place the sterile normal saline bag into the pressure bag. Inflate the pressure bag to between 150 and 300 mm Hg. Purge the system of air until fluid is seen coming from the openings on the distal tip of the device (Figure 2). Clamp off tubing set until ready for use.

Figure 2

7. Attach the pen electrical cable to the Cardioblate generator connector per instructions outlined in "Connections to External Devices."
8. Prior to tissue ablation:

- a. Open the IV line and make certain that saline is flowing freely out of the openings in the distal tip of the device. Continuous flow of saline is necessary during ablation to cool the outer tissue. Cooling of the outer tissue prevents charring and tissue desiccation and allows radiofrequency energy to be delivered deeper into the tissue to more effectively create transmural lesions.
 - b. Set the radiofrequency ablation generator to the POWER mode and set the power between 20-30 watts maximum output. Ideal power settings should be determined by the surgeon during tissue ablation.
 - c. Position the device electrical cable to avoid contact with the patient or other leads.
9. The distal tip of the device is dragged gently across the targeted intracardiac tissue under direct visualization by the surgeon. The current is activated in the device using the foot switch provided with the Cardioblate generator. The current flows between the ablation tip and the large indifferent electrode on the patient's body.
10. The lesion becomes visible as a white coloration on the tissue. The tip is moved in an oscillating motion at a recommended speed of 1 cm/sec. The surgeon will determine depth, width, and location of the ablation.

Table 1. Anticipated Lesion Depth at Various Power Levels

95% Confidence interval for 10 second oscillation, 1 cm	95% Confidence interval for 20 second oscillation, 1 cm	Power (W) @ 150 - 300 mm Hg pressure on saline irrigation
2.3 - 3.5 mm	3.5 - 4.4 mm	20
2.9 - 3.7 mm	4.7 - 5.7 mm	30

11. Do not dwell in a target area too long; sudden tissue "popping" and/or perforation may occur.
12. Upon completion of the surgical procedure, disconnect the device and indifferent electrode from the Cardioblate generator and discard.

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Disclaimer of Warranty

ALTHOUGH THE CARIOBLATE XL SURGICAL ABLATION PEN, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law, by a court

of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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Disclaimer of Warranty

ALTHOUGH THE CARIOBLATE XL SURGICAL ABLATION PEN, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

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© 2018 Medtronic
Trademarks may be registered and are the property of their respective owners.