« 22 » September 2021

To be submitted to the Federal Service for Surveillance in Healthcare (Roszdravnadzor) in Russian Federation.

We, Olympus Winter & Ibe GmbH, located at Kuehnstr. 61, 22045 Hamburg, Germany, as the legal manufacturer of the Medical device "Electrosurgical generator ESG-410 with accessories", hereby state that attached documents is true:

1. Electrosurgical generator ESG-410 Instruction for use including operating instructions and a list of contact details for the client
2. Quick reference guide "Sockets and screens"
3. Quick reference guide "Footswitches, neutral electrodes and procedures"
4. Quick reference guide "Mode lookup table"
5. Maintenance manual for regular safety checks
6. Single-pedal switch WA91303W and Two-pedal switch WA91301W instruction manual and contact information for the client
7. Two-pedal switch, wireless WA91302W instruction manual and contact information for the client
8. Electrosurgical generator ESG-410, Single-pedal switch WA91303W, Two-pedal switch WA91301W, Two-pedal switch, wireless WA91302W addendum to instructions for use

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Dr. Christian Meyer
Managing Director

OLYMPUS SURGICAL TECHNOLOGIES EUROPE

Olympus Winter & Ibe GmbH, Kuehnstraße 61, 22045 Hamburg, Germany, P.O. Box 70 17 09, 22017 Hamburg, Germany

Phone: +49 40 669 66-0, Fax: +49 40 669 66-2109, www.olympus-oste.eu

Managing Directors: Dr. André Roggan (Executive Managing Director), Kazutaka Eguchi, Dr. Christian Meyer, Tomohisa Sakurai, Akihiro Taguchi, Carl Constantin Zangemeister, Reinhard Zentner


Trade Registration Court: Amtsgericht Hamburg HRB 16 328

I, the undersigned Dr Matthias Kleiser, Civil Law Notary of this Free and Hanseatic City of Hamburg, Federal Republic of Germany, hereby certify and attest that

Mr. Dr. Christian Martin Meyer
date of birth: 17th day of June 1977
residing: Kuehnstr. 61, 22045 Hamburg
personally known to me, the Notary,
personally signed before me.

IN WITNESS WHEREOF I have hereonto set my hand and affixed my seal of office.

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
CE 0197

WA91307C

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Contents

1 Using this document.....	7
1.1 User instructions	7
1.2 Conventions throughout this document.....	7
1.3 Safety-related information in this document	7
2 Safety information.....	9
2.1 Intended purpose	9
2.2 Contact.....	10
2.3 General safety notes	10
3 Product description	19
3.1 Scope of delivery	19
3.2 Product functions	19
3.3 Operating elements on the hardware	21
3.4 Operating elements on the touchscreen	22
3.5 Popups and messages.....	25
3.6 Symbols and icons	26
3.7 Warranty	29
4 Installation	30
4.1 Safety notes for installation	30
4.2 General inspection	30
4.3 Placement of the electrosurgical generator.....	30
4.4 Connection of electrosurgical generator to the mains electricity	31
4.5 Connection of foot switches	32
4.6 Connection of peripheral equipment	32
4.7 Connection for automatic smoke evacuation	32
5 System settings.....	36
5.1 Operation of the touchscreen.....	36
5.2 Settings – Setup.....	36
5.3 Volume control of activation tones	37
5.4 Brightness control	38
5.5 Alarm tones	38
5.6 Touch tone	39
5.7 Languages	39
5.8 Date, time and format.....	40
5.9 Monopolar modes	40
5.10 Startup settings	41
5.11 Procedures.....	41
5.12 Assigning foot switches.....	42
6 Safe use of neutral electrodes	44
6.1 Split type and non-split type neutral electrodes and CQM	44
6.2 Restriction to only split type neutral electrodes	45
6.3 Safety feature for baby neutral electrodes	46
6.4 Conditions of the CQM indicator	48
6.5 Selecting an appropriate neutral electrode	49
6.6 Connection of a neutral electrode to the electrosurgical generator.....	49
6.7 CQM indicator verification	51

7 Before use.....	52
7.1 Safety notes for before use	52
7.2 Inspection before use.....	53
7.3 Overview of modes	55
7.4 Description of modes	55
7.5 Tissue effects depending on the power level and effect level.....	62
8 Use.....	64
8.1 Safety notes for use	64
8.2 Description of the output sockets	65
8.3 Connection of instruments	67
8.4 Selection of appropriate output settings.....	67
8.5 Automatic instrument recognition.....	69
8.6 Selecting a user-defined procedure	72
8.7 Activation of HF output.....	72
8.8 Activating.....	73
8.9 Alarm system	77
8.10 Troubleshooting.....	78
9 After use	83
10 Reprocessing of the ESG-410.....	84
10.1 General information for reprocessing.....	84
10.2 Cleaning.....	84
10.3 Disinfection.....	84
10.4 Reprocessing of other HF equipment	85
11 Maintenance	86
11.1 Regular safety checks.....	86
11.2 Maintenance settings	87
12 Repair and shipment.....	92
12.1 Repair	92
12.2 Shipment.....	92
13 Storage and disposal.....	93
13.1 Storage.....	93
13.2 Disposal	93
14 Ordering information	94
15 Compatible equipment	95
15.1 System chart	95
15.2 Compatible neutral electrodes	96
15.3 Compatible foot switches	97
16 Technical data	98
16.1 Specifications for the ESG-410.....	98
16.2 Ambient conditions for operation, storage and transport	98
16.3 Plug requirements	99
16.4 Resistance range for CQM.....	99
16.5 Specifications for cables	99
16.6 Specifications for the USB flash drive	99
17 Additional information.....	100
17.1 Alarm system	100

17.2 Tone information	101
17.3 Mode characteristics according to IEC 60601-2-2	104
17.4 Characteristics of High Power Cut Support (HPCS)	105
17.5 Output characteristic diagrams	105
18 Electromagnetic compatibility	132
19 Glossary.....	135
20 Index.....	139

1 Using this document

This chapter provides general information on how to use these instructions for use.

1.1 User instructions

Before use, thoroughly read these instructions for use and the instructions for use of all other products that will be used during the procedure.

If the required instructions for use are missing, immediately contact an Olympus representative.

Keep the instructions for use in a safe, accessible location.

1.2 Conventions throughout this document

Element	Explanation
	This is the safety alert symbol. It is used to alert the user to potential physical injury hazards. Observe all safety messages that follow this symbol to avoid possible injury.
	This symbol indicates additional helpful information.
1. 2. 3.	A numeration indicates a sequence of actions.
•	Bullet points indicate individual actions.
►	Arrows indicate the results of actions.
— — —	Dashes indicate the listing of data, options or objects.
1) 2) 3)	Numbers with right parenthesis name elements in illustrations.
[...]	Bracketed terms refer to elements in the graphical user interface or keys. Elements in the graphical user interface can be: - buttons - menu items - dialog elements

1.3 Safety-related information in this document

For safe and efficient use of this product it is important to read the complete instructions for use. The safety related information of these instructions for use is not only described in this chapter. All information and instructions in this document are to ensure a safe and efficient use.

Safety information are described in 3 categories:

- General safety notes
- Specific safety notes
- Safety messages

1.3.1 General safety notes

General safety notes are located at the beginning of the instructions for use. They address safety issues that apply to the general use and handling of the product.

1.3.2 Specific safety notes

Specific safety notes are located at the beginning of a chapter or section. They address safety issues that are specific to the information and actions described in the chapter or section.

1.3.3 Safety messages

Safety messages focus on unexpected incidents that can occur during the intended use. They are placed in the context of such risks and actions. The safety messages are highlighted by a signal word.

WARNING

Indicates a potentially hazardous situation which, if not avoided, can result in death or serious injury.

CAUTION

Indicates a potentially hazardous situation which, if not avoided, can result in minor or moderate injury.

NOTICE

Indicates a property damage message.

2 Safety information

This chapter provides general safety information for using the product. The safety information includes:

- Intended purpose
- Manufacturer information
- General safety notes

2.1 Intended purpose

2.1.1 Intended use and indications

The electrosurgical generator, in conjunction with electrosurgical accessories and ancillary equipment, is intended for cutting and coagulation of tissue in the following medical fields:

- Open surgery
- Laparoscopic surgery
- Endoscopic surgery

Only for use by a qualified physician in an adequate medical environment.

2.1.2 Contraindications

Absolute contraindications

There are no absolute contraindications.

Relative contraindications

At the discretion of the physician, electrosurgical procedures can be contraindicated under the following circumstances:

- If tissue coagulation and cutting can have a negative effect on the state of the patient.
- For patients with implanted electronic devices, e.g. cardiac pacemaker, cardioverter defibrillators.
- For patients with a weakened immune system.
- For patients with blood coagulation disorders.

2.1.3 Patient target group

Not dedicated to a specific patient population. Patient target groups are defined by the individual applicator used with the device.

2.1.4 Intended users

Medical use

Only for use by a qualified physician or medical personnel under the supervision of a physician.

Only for use in an adequate medical environment.

These instructions for use do not explain or discuss clinical procedures.

Reprocessing

Reprocessing of reusable products may only be performed by qualified hygiene personnel.

Repair

Repair of the product must only be performed by trained qualified servicing personnel that has been authorized by Olympus. Otherwise, Olympus cannot be held responsible for the safety and performance of the product.

Be aware that also fuses must only be changed by an authorized service center.

Installation and setup

The electrosurgical generator and its peripheral equipment must be properly installed and set up by qualified personnel that has been authorized by Olympus.

Regular safety check

The regular safety check must only be performed by a qualified electrician with sufficient experience in maintaining medical electrical devices.

The regular safety check must only be performed in a technical laboratory with appropriate equipment.

Software update

Software updates must only be performed by qualified personnel.

2.1.5 Clinical benefits

Clinical benefits depend on the applied part used with the device.

2.2 Contact

2.2.1 Manufacturer

Olympus Winter & Ibe GmbH
Kuehnstr. 61
22045 Hamburg
Germany

Country-specific information

Refer to the contact information for customers included in the scope of delivery for country-specific information.

2.2.2 Incident reporting

If a serious incident occurs with the device, report it to the manufacturer and/or the relevant national authority according to national legislation.

2.3 General safety notes

The following safety information address points that apply to the general handling and use of the product. This information is to be supplemented by the information given in each chapter of this document.

Read all instructions for use thoroughly

Before use, read this instructions for use and the instructions for use for all other equipment which will be used during the procedure.

An insufficient understanding of the safety information in these instructions for use can result in death, serious injury, or equipment damage.

Liquids and splashes

If liquids get into the electrosurgical generator, there is a risk of electric shock and damage to the product.

Environmental conditions

Inappropriate environmental conditions during use, storage or transport can damage the electrosurgical generator and its accessories.

Vibration and shocks

Vibrations and shocks during use, storage or transport can damage the electrosurgical generator and its accessories.

Appropriate power supply and grounding

If the electrosurgical generator is improperly connected to the mains electricity there is a risk of electric shock.

Potential equalization

The purpose of additional equipotential bonding is to equalize potentials between different metal parts that can be touched simultaneously in the medical environment. Equipotential bonding increases the electrical safety.

Compatible equipment

Using incompatible equipment can lead to excessive or incorrect HF output. This can result in injury to the patient and the user as well as damage to the products.

Secure connection of the plugs

All plugs must be securely inserted to the corresponding sockets. Be aware that the electrosurgical generator immediately turns off if the mains plug is accidentally disconnected. During a procedure, this can result in complications for the patient.

If the foot switch plug is not connected securely, then the HF output may not be activated. There is the risk that the HF instrument cuts tissue mechanically. This can result in bleeding and perforation of the tissue.

Electromagnetic compatibility

During activation of the electrosurgical generator, high frequency signals or spark discharge noise can disturb neighboring electrical equipment. Malfunction of the devices can occur, e.g. the monitor of endoscopic imaging equipment can freeze or black out. The precautionary measures are described in the chapter "18 Electromagnetic compatibility" on page 132. Additionally, the instructions regarding installation and connection in the following chapters must be observed.

Damaged equipment and accessories

The use of damaged equipment, of equipment with improper functioning or the reuse of single use products can cause electric shock, mechanical injury, infection and thermal injury to the patient and the user.

Procedural hazards and complications

If a malfunction of the electrosurgical generator occurs during the procedure, there is the risk that the procedure is interrupted, delayed or cannot be completed. A prolonged use or an energy dosage more than necessary can compromise the patient safety.

Disconnection of cables

Tugging and pulling on the cables can break the electric circuit and can damage the equipment.

Maintenance

To ensure electrical safety and to maintain device performance, Olympus recommends to perform a safety check annually or according to national regulations and guidelines on all medical electrical equipment.

Responsibility of the user

In electrosurgery, electric current is applied to a patient to modify, destroy, dissect or remove tissue. This method bears certain risks for the patient, the operator and the assisting medical personnel. It is important that all users are aware of such risks.

The information in the following sections summarize important basics and rules of electrosurgery to remind the operator and the medical personnel of the associated risks. This information is not meant to replace the fact that the all users need to have a profound understanding of electrosurgical procedures.

2.3.1 Fire and explosion

Sparks that are generated in most electrosurgical procedures can ignite surroundings, substances and materials. Fire and explosion can occur.

Flammable anesthetics and fire supporting gases

Especially consider:

- Flammable anesthetics or oxidizing gases, e.g. nitrous oxide (N₂O) or oxygen
- Volatile solvents, e.g. ether or alcohol

Flammable agents for cleaning and disinfection

Non-flammable agents should be used for cleaning and disinfection wherever possible. If flammable agents are used, they should be allowed to evaporate before the application of HF surgery. There is a risk of pooling of flammable solutions under the patient or in body depressions such as the umbilicus, and in body cavities such as the vagina. Any fluid pooled in these areas should be mopped up before HF surgical equipment is used.

Flammable gases in the gastrointestinal tract/endogenous gases

When performing electrosurgery in the gastrointestinal tract, flammable gases can be present in the intestines.

Attention should be called to the danger of ignition of any endogenous gases.

Flammable materials

An active electrosurgical instrument or an instrument that is still hot from use can cause fire when in contact with flammable materials, e.g. gauze, surgical drape or body hair.

2.3.2 Electrosurgical smoke

The thermal destruction of tissue creates a smoke by-product. This smoke can contain toxic gases and vapors such as benzene, hydrogen cyanide, formaldehyde, bioaerosols, dead and live cellular material (including blood fragments) and viruses.

- At high concentrations the smoke causes ocular and upper respiratory tract irritation.
- The smoke creates visual problems for the surgeon.
- The smoke has been shown to have mutagenic potential.

- The smoke has unpleasant odors.

Olympus recommends to take adequate measures during the procedure to protect the surgical staff from inhaling unfiltered smoke, e.g. by using smoke evacuators, sufficient ventilation in the operating room or wearing surgical masks.

2.3.3 Effects of electric current on biological tissue

Electrolytic effect – undesired

When direct current (DC) is applied to biological tissue, spatial separation of positive and negative ions occurs. This produces acids and bases at the electrode poles and results in acid burns to the tissue. Therefore, only alternating current (AC) is used for electrosurgical procedures.

Faradic effect or neuromuscular stimulation – undesired

Alternating current (AC) that is provided by the mains electricity has a frequency of approximately 50 to 60 Hz. Electric current at such a frequency has undesired effects to the human body like pain, nerve stimulation, muscle stimulation or cardiac arrest. As the human body is less sensitive to current with high frequencies, electrosurgical generators transform the low frequency AC provided by the mains electricity to high frequencies of above 300 kHz, i.e. high frequency current or HF current.

However, due to high current densities, low frequency currents can still arise when using an electrosurgical generator. Reasons for this phenomenon are, e.g.:

- Spark discharge to the tissue or to other metallic objects
- High HF output settings
- The use of modes with HPCS (High Power Cut Support), see the section "7.4 Description of modes" on page 55.

Therefore, it's important to be aware that neuromuscular stimulation - which can cause pain and discomfort to patients without sedation, pain medication or general anesthesia - can occur during electrosurgery.

Thermal effect – desired

When the applied HF current flows through the tissue and needs to overcome its electric resistance, heat is generated within the tissue. In other words, it is not the metal part of the electrosurgical instrument that heats the tissue, it is the electric current that flows within the tissue. This thermal effect of the electric current on tissue is used for cutting and coagulating.

Depending on the specific resistance of different tissue types, the thermal effect varies. The higher the tissue resistance is, the less amount of current can flow through the tissue and the smaller the desired thermal effect is.

Basically, tissue that has a high water content has a low resistance and therefore the thermal effect is high. Some examples of specific resistances in different tissue types:

- Muscle: 160 to 250 Ω/cm
- Blood: 160 to 300 Ω/cm
- Liver: 200 to 380 Ω/cm
- Brain: 670 to 700 Ω/cm
- Fat: 1,600 to 3,300 Ω/cm
- Skin: 11,000 Ω/cm

Compared to metals, e.g.:

- Copper: 0.0000017 Ω/cm

- Silver: 0.0000016 Ω/cm

2.3.4 Electrosurgical effects – coagulation, cutting and carbonization


Coagulation

A temperature of 60 to 70 °C in the area around the active electrode gradually heats the tissue cells until their intracellular fluid escapes through the cell membranes. The results are:

- Protein denaturation
- Shrinking of the cells
- Shrinking of vessels
- Clotting of blood
- Hemostasis

Cutting

When the generated heat causes the intracellular fluid to rapidly evaporate, then this results in an explosive burst of the cell membrane. To achieve this, a very high voltage of approximately 200 V_p is required. Such a high voltage generates sparks between the tissue and the active electrode. These sparks concentrate the HF current to specific points of the tissue that heats up to more than 100 °C. The fast repetition of this effect around the active electrode determines the cut.

Principle of an electrosurgical cut

It is important to understand that the sparks are only able to ignite repetitively if the active electrode is not in direct contact with the tissue.

The desirable side effect of electrosurgical cutting is the hemostasis of neighboring tissue cells, as part of the heat is transferred to surrounding tissue where coagulation takes place.

Risk of mechanical injuries

To avoid the risk of mechanical cuts by the metallic tip, blade or snare of an HF instrument it is important to keep in mind that it is the HF current that performs the desired cut. If too much mechanical stress is applied to the target tissue, e.g. by a loop that is fastened to strongly around the target tissue, then bleeding or perforation can occur.

Carbonization – particularly undesirable

Carbonization is the overheating of desiccated tissue. The remains of the tissue cells are reduced to carbon. Carbonization of tissue takes place from approximately 200 °C upwards. It is important to understand that tissue carbonization must be avoided as it impairs the wound healing.

2.3.5 The electric circuit**Three important rules of electricity**


When working with electric energy the following 3 basic rules must be respected:

- Electric current needs a complete circuit to flow.
- Electric current is 'lazy' and prefers the path of least resistance.
- Electric current seeks ground.

The circuit of electrosurgical HF current

In electrosurgery the complete electric circuit consists of:

- The electrosurgical generator
- A monopolar or bipolar electrosurgical instrument for cutting or coagulation
- The patient's target tissue
- For monopolar procedures: a neutral electrode
- The connecting cables

2.3.6 Monopolar electrosurgery

In monopolar electrosurgery the HF current passes from the active electrosurgical instrument through the patient to the neutral electrode and back to the electrosurgical generator to complete the electric circuit.

High current densities and thermal effect at the instrument's tip

The electrosurgical instrument usually has a pointed shape with a small surface. Therefore, a high current density accumulates at the instrument's tip which generates enough heat for cutting or coagulating the target tissue.

Low current densities and thermal effect at the neutral electrode

The neutral electrode has a large surface so that the current is spread over a large area. This reduces the current density and the heat. As the electric current always takes the path of the lowest resistance, the neutral electrode ensures that the current returns to the electrosurgical generator without causing overheating of tissue outside the surgical site.

The entire surface of the neutral electrode must be reliably attached to the prepared patient's skin, usually on the thigh.


Neutral electrodes

Wrong application, wrong size or malfunction of the neutral electrodes are the main causes for endogenous burns.

- Observe the information in the chapter "6 Safe use of neutral electrodes" on page 44.

2.3.7 Bipolar electrosurgery

For surgical procedures where the HF current could flow through parts of the body that have a relatively small cross sectional area, the use of bipolar techniques can be desirable to avoid unwanted tissue damage.

In bipolar electrosurgery no neutral electrode is required. Bipolar electrosurgical instruments are designed in a way that the current passes only between the 2 electrodes that are implemented within the same device, e.g. between the jaws of a bipolar forceps. In bipolar electrosurgery the current does not pass large areas of the patient's body but only through the tissue that is between the two electrodes within the electrosurgical instrument.

The bipolar technique makes electrosurgery much safer, because:

- The current exposition to the patient is limited.
- There is no risk of burns caused by the neutral electrode. Refer to the section "6 Safe use of neutral electrodes" on page 44.
- The risk of leakage currents is reduced. Refer to the section "2.3.9 Unintended current flows" on page 17.

2.3.8 Implants

Electronic implants like cardioverter defibrillators, cardiac pacemakers or cochlear implants are very sensitive to electric energy. Therefore, electrosurgical procedures can interfere with and can damage such devices. Also other electrically conductive implants like metal implants present a hazard due to concentration or re-direction of HF currents. Qualified advice should be obtained before the electrosurgical procedure.

2.3.9 Unintended current flows

Physiological monitoring equipment

When using physiological monitoring during an electrosurgical procedure, current can flow to the monitoring electrodes. To prevent the risk of high current densities, the monitoring electrodes need to be placed as far as possible from the surgical electrodes (instrument, neutral electrode). Especially monitoring electrodes in form of needles bear a high risk of burns at the site where the electrodes are attached. Therefore it is recommended to not use needle monitoring electrodes. In any case it is recommended to use monitoring equipment with HF current limiting measures.

Conductive fluids within the patient

Conductive fluids within the patient, like blood or saline, that are in direct contact with or in close proximity to the active instrument can transmit electric current and heat to non-target tissue. This can cause unintended burns to the patient.


Earth leakage currents

Unfortunately, alternating HF current techniques always correspond with a certain level of earth leakage currents. Earth leakage currents are current flows that pass another pathway than the initially intended electric circuit. A certain percentage of the current could pass via the patient and the operating table or the surgeon and then via the ground back to the electrosurgical generator. Earth leakage currents bear the risk of thermal tissue damages.

Olympus electrosurgical generators constantly monitor earth leakage currents and give an acoustic and visual feedback if the earth leakage current exceeds a certain level. These feedbacks make the physician aware that there is an increased risk of burns for the patient.

The user must not touch the patient and any contacts in the sockets simultaneously.


Capacitive coupling

Another phenomenon of electricity is that alternating HF current can pass from a conductive material through an insulator into another nearby conductive material without the aid of wires or cables. This can occur despite intact insulation surrounding the conductors and is called capacitive coupling. So it is possible e.g. that the HF current passes via the insulation of an active electrosurgical instrument to a metallic trocar which is in contact with some tissue in the bowel. If the contact area is small so that the current density is high, then this can result in a burn to the bowel tissue.

Capacitively coupled currents mainly occur in monopolar electrosurgery. They are also called leakage currents or stray currents.


Direct coupling

Direct coupling is the current flow directly from one conductor to another conductor. For example, when an active electrosurgical instrument like a cutting electrode is in contact with another metallic instrument, like a grasping forceps. Then current flows from the active cutting electrode into the grasping forceps and the grasping forceps can cause burns to non-target tissue.

For endoscopic surgery it is important to understand that the electrosurgical instrument must only be activated when the conductive (non-insulated) part of the active electrode is in full view.

Insulation failures

When the insulation of an electrosurgical instrument or of an HF cable is damaged the electric current can pass another way than the intended way. If the insulation is damaged, another metal instrument can be activated or non-target tissue inadvertently can be burned.

3 Product description

This chapter provides a general overview of the product and its functions. This includes the following information:

- Scope of delivery
- Product functions
- Operating elements
- Explanation of symbols and icons
- Warranty information

3.1 Scope of delivery

Before use, check that all items listed below are available.

Contact an Olympus representative or an authorized service center if any items are missing or damaged.

- Electrosurgical generator ESG-410
- Instructions for use, including contact information sheet for customers
- Quick Reference Guide "Sockets and screens"
- Quick Reference Guide "Foot switches, neutral electrodes and procedures"
- Quick Reference Guide "Mode lookup table"
- Maintenance Manual for the regular safety check

3.2 Product functions

3.2.1 Operating principles

The electrosurgical generator ESG-410 is designed to be used for the following electrosurgical applications:

- Bipolar electrosurgical cutting and coagulation.
- Monopolar electrosurgical cutting and coagulation, which requires a neutral electrode.

The several cutting and coagulation modes are designed for procedures in open, laparoscopic and endoscopic surgery. All settings of the device are selected via the touchscreen on the front panel.

The ESG-410 is activated either via a compatible foot switch or via a hand switch at the connected instrument.

The ESG-410 is connected to the local power supply.

Details on the hardware and software, on the correct handling and use of the device are described in the following chapters.

3.2.2 Features

The electrosurgical generator is equipped with the following features:

- Contact Quality Monitor (CQM)
- High Power Cut Support (HPCS)
- Fast Spark Monitor (FSM)

CQM – Contact Quality Monitoring for the neutral electrode

When using split type neutral electrodes for monopolar electrosurgery, the ESG-410 is able to detect unintended detachment of the neutral electrode from the patient. The indicator for the contact quality monitor (CQM) is green while the contact between the split type neutral electrode and the skin of the patient is within an acceptable resistance range. If the contact between the split type neutral electrode and the patient's skin is insufficient, then an alarm tone sounds, a warning message is displayed and the CQM indicator is red.

Using non-split type neutral electrodes is not as safe as using split type ones. CQM is not able to detect any detachment of non-split type neutral electrodes from the patient.

For detailed information on the safe and correct use of neutral electrodes refer to the chapter "6 Safe use of neutral electrodes" on page 44.

HPCS – High Power Cut Support


This feature optimizes the start of the cutting procedure for certain cutting modes. By applying high power to the tissue, HPCS supports immediate spark ignition and reduces the risk of mechanical cutting.

FSM – Fast Spark Monitor

When using cutting modes this feature ensures smooth and reproducible cutting although the tissue characteristics are varying, e.g. in muscle and fat.

3.3 Operating elements on the hardware


3.3.1 Front panel

- 1) Output socket UNIVERSAL 1*
To connect an Olympus HF instrument with instrument recognition.
- 2) Output socket UNIVERSAL 2*
To connect an Olympus HF instrument with instrument recognition.
- 3) Output socket BIPOLAR*
To connect a bipolar HF instrument.
- 4) Output socket MONOPOLAR 1*
To connect a monopolar HF instrument.
- 5) Touchscreen
To display the status of connected accessories.
To show and modify settings.
- 6) Output socket MONOPOLAR 2*
To connect a monopolar HF instrument.
- 7) Neutral electrode socket*
To connect a neutral electrode.
- 8) Power switch
To switch the electrosurgical generator on and off.

*Applied part according to standard IEC 60601-1.

3.3.2 Rear panel


- 1) Volume control
To adjust the volume of the activation tone.
- 2) Ventilation slots
- 3) Potential equalization connector
To increase electrical safety by potential equalization
- 4) Fuse holder
Compartment for the replaceable fuses.
- 5) Power socket
To connect the ESG-410 to the mains electricity.
- 6) Foot switch sockets
To connect a compatible foot switch.
- 7) LINK-OUT socket
To connect a communication cable for peripheral equipment with a 14-pin plug.
- 8) LINK-IN socket
To connect a communication cable for peripheral equipment with a 26-pin plug.
- 9) USB socket
To connect a USB flash drive.

3.4 Operating elements on the touchscreen

The control buttons on the touchscreen are activated by tapping the corresponding part on the screen with the finger tip. The home screen displays the selected settings for the different output sockets. When tapping the pane for one of the output sockets the display changes to the set screen.

3.4.1 Home screen

Each of the 6 panes on the home screen corresponds to one of the 6 output sockets of the ESG-410.

UNIVERSAL panes accessible after connecting instrument

The 2 Universal panes on the left side represent the 2 Olympus specific output sockets UNIVERSAL 1, UNIVERSAL 2. These 2 panes are accessible only when a correspondent Olympus instrument is connected to the output socket.

- 1) Selected user-defined procedure
- 2) Toggle key
To switch the output socket on and off.
- 3) Selected power setting
- 4) Selected effect setting
- 5) Socket pane (switched off)
- 6) Selected coagulation mode
- 7) [Settings] button
To enter the setup and maintenance screen.
- 8) [Procedure] button
To select or save a user-defined procedure.
- 9) [Foot switch] button
To assign a foot switch to an output socket.
- 10) [COM] button
- 11) [Info] button

10) CQM info button

The neutral electrode indicator displays the current state of the CQM.

The information screen explains the different states of the CQM.


The neutral electrode setup provides the setting for baby neutral electrodes.

11) Selected cutting mode

12) Foot switch indicator

3.4.2 Set screen

On the set screen the required modes and power levels are selected.

- 1) Foot switch indicator
- 2) Selected user-defined procedure
- 3) Coagulation mode button
Displays the selected coagulation mode.
Tap to change the mode.
- 4) Selected power level
Tap + or - to change the power level.
- 5) Selected effect level
Tap the arrow to change the effect level.
- 6) Mode information button
Tap to receive information on the modes.
- 7) Cutting mode button
Displays the selected cutting mode.
Tap to change the mode.

3.5 Popups and messages

The ESG-410 displays 3 different categories of messages:

- Error messages
- Info messages
- Message boxes

Error messages


These messages are displayed when the system detects a safety relevant issue, e.g. if the connections between the ESG-410 and a peripheral equipment like a videoprocessor are not correct. Error messages provide instructions (remedial actions) to clear the alarm.

When an error message pops up, an alarm tone sounds. More details about error messages are described in the section "8.9 Alarm system" on page 77.

Info messages

These kind of messages assist the user if something has not been done correctly, so that the system cannot work properly. When an info message pops up, an info tone sounds. The message can be cleared by tapping the [OK] button.


If, for example, the user has missed to assign the foot switch, the following info message can pop up:

- 1) Information symbol and code
- 2) Message title
- 3) Remedial actions
- 4) [OK] button to clear the message

Message boxes

Message boxes provide helpful information to assist the user during a certain operation. For example, the following message box pops up when a new foot switch is connected to the ESG-410:

- 1) Message title
- 2) Remedial actions
- 3) Confirmation button to clear the message


When a message box pops up, an info tone sounds. The message can be cleared by tapping the confirmation button, e.g. [OK] or [Cancel].

3.6 Symbols and icons

Symbols

The following table gives an explanation for each symbol used on the product and on the packaging of the product.

Symbol	Explanation
SN	Serial number
QTY	Quantity
REF	Catalog number

Symbol	Explanation
Rx only	Federal (USA) law restricts this device to sale by or on the order of a physician
	Medical device
	Refer to instructions for use
	Caution, consult accompanying documents
	Defibrillation-proof type CF applied part
	Manufacturer
	CE marking indicating that the device is in conformity with the applicable requirements set out in applicable European Union harmonization legislation
	Unique device identifier
	Transport conditions
	Storage conditions
	Temperature limit
	Humidity limit
	Keep dry
	General symbol for recovery/recyclable
	Green Dot symbol for dual recycling system
	In accordance with European Directive on Waste Electrical and Electronic Equipment (WEEE), this symbol indicates that the product must not be disposed of as unsorted municipal waste, but should be collected separately
	Fragile, handle with care
	This way up
	Stacking limit by number
	HF isolated patient circuit; symbol for patient circuit where there are no components installed to provide a low-impedance path to earth for HF currents
	Volume control

Symbol	Explanation
	USB socket
	Foot switch
	Certification mark of the TUV Rheinland of North America, Inc.
	Fuse
	Equipotentiality terminal
	Non-ionizing radiation; indicates that the device includes RF transmitters or that the device intentionally applies RF electromagnetic energy for diagnosis and treatment

Icons

The following table gives an explanation for each icon that can be displayed on the touchscreen of the electrosurgical generator.

Icon	Explanation
	AUDIO ON An alarm tone sounds for errors with low or medium priority.
	AUDIO OFF No alarm tone sounds for errors with low or medium priority.
	Autostart ON Activation starts automatically with tissue contact.
	Autostart OFF Activation does not start automatically with tissue contact.
	Single-pedal foot switch, wired
	Double-pedal foot switch, wired
	Double-pedal foot switch, wireless
	There is no neutral electrode connected.
	CQM detects insufficient contact between the neutral electrode and the patient's skin.
	CQM detects that the contact between the neutral electrode and the patient's skin is fine.
	There is a non-split type neutral electrode connected. CQM cannot detect insufficient contact between the neutral electrode and the patient's skin.
	The safety setting for baby neutral electrodes is switched on. There is no neutral electrode connected.
	The safety setting for baby neutral electrodes is switched on. CQM detects insufficient contact between the neutral electrode and the baby's skin.

Icon	Explanation
	The safety setting for baby neutral electrodes is switched on. CQM detects that the contact between the neutral electrode and the baby's skin is fine.
	The safety setting for baby neutral electrodes is switched on. There is a non-split type neutral electrode connected. CQM cannot detect insufficient contact between the neutral electrode and the baby's skin.
	Maintenance reminder The due date for the next regular safety check is reached.

3.7 Warranty

Any warranty claims towards Olympus are forfeited if the user or unauthorized persons attempt repair or modification of the product. No warranty is provided for any damage due to misuse of the product.

4 Installation

This chapter describes all necessary steps to install the product at the installation site before its first use.

4.1 Safety notes for installation

An improper installation of the electrosurgical generator and its equipment compromises the safety. To avoid injuries and equipment damage observe the instructions for installation in the following sections.

4.2 General inspection

Before installing and connecting, inspect the electrosurgical generator and its equipment.

- Check that the electrosurgical generator, its equipment and all accessories do not have:
 - any dents, cracks, kinks, or deformations
 - any deep scratches
 - any corrosion
 - any missing or loose parts
- Check all markings on the electrosurgical generator, its equipment and all accessories for clear visibility.

4.3 Placement of the electrosurgical generator

General requirements for correct placement

To prevent malfunction or damage of the electrosurgical generator and its accessories, consider the following circumstances when placing the electrosurgical generator:

- Place the electrosurgical generator on a stable, level surface.
- Put the electrosurgical generator upright onto its feet. Do not install the electrosurgical generator on its side or upside down.
- If the electrosurgical generator is placed on a trolley, the trolley must be of adequate strength and size to hold the electrosurgical generator securely.
- The compatible trolley for the ESG-410 is the Energy Cart TC-E410.
- Make sure that there is sufficient space at the back and at the sides of the electrosurgical generator. This allows air circulation at the ventilation slots.
- Only use the electrosurgical generator in compliance with the environmental conditions as specified in the chapter "16 Technical data" on page 98.
- When lifting up, grip the electrosurgical generator at its sides with both hands. Do not hold the electrosurgical generator at its fuse holder.

Electromagnetic considerations

- Make sure that the electrosurgical generator is not used adjacent to or stacked with equipment that is not part of this electrosurgical generator or system.
- Thoroughly confirm the compatibility of all equipment.
- Do not use the electrosurgical generator in conjunction with:
 - Electrical equipment for which the safety against leakage current is not confirmed.
 - Electrosurgical equipment for which the safety in combined use is not confirmed.
- Do not loop cables.
- Do not bundle cables together with cables belonging to other medical equipment.
- Follow the instructions regarding electromagnetic ambient conditions given in the chapter "18 Electromagnetic compatibility" on page 132.

4.4 Connection of electrosurgical generator to the mains electricity

WARNING **Electric shock**

If an unknowingly damaged electrosurgical generator is connected inadequately to the mains electricity, there is the risk that voltage is applied to the housing by internal electrical components. This can result in electric shock to the user.


Only connect the electrosurgical generator to a supply mains with protective earth (grounding).

Before connecting

- Make sure that the fixed socket outlet is grounded. Only connect the electrosurgical generator to a grounded fixed socket outlet.
- Make sure that the fixed socket outlet meets the power requirements. Refer to the section "16 Technical data" on page 98 or to the identification plate on the back panel of the electrosurgical generator.
- Only use power cables as specified in the section "16 Technical data" on page 98.
- Only connect the electrosurgical generator if the power specifications of the fixed socket outlet match the specifications.
- If a multiple socket outlet is used, consider the following:
 - The multiple socket outlet must be equipped with an insulating transformer of protection class I, conforming to IEC 60601-1.
 - Observe the maximum permitted current or power loading of the multiple socket outlet and the insulating transformer.
- Make sure that the mains plug is securely connected to the socket outlet.
- Do not place portable socket outlets on the floor to avoid trip hazards which can unplug the power cable.
- Make sure that the fuses are sufficient for the power requirements of the electrosurgical generator and the additional equipment.
- Make sure to have a compatible power cable available.

Connecting

1. Make sure that the power switch of the electrosurgical generator is not depressed.

2. Connect one plug of the power cable to the power socket on the rear panel of the electrosurgical generator.
3. Connect the other plug of the power cable to the socket outlet of the mains electricity.
4. Connect the potential equalization connector to an equipotential bonding lead.
5. Position the electrosurgical generator in a way that the power cable can be disconnected from the mains electricity.

4.5 Connection of foot switches

Up to 3 foot switches can be connected to the rear panel of the electrosurgical generator. However, for safety reasons it is not possible to use foot switches of the same type simultaneously, e.g. 2 single-pedal foot switches or 2 wireless double-pedal foot switches. Therefore, ensure that only different types of foot switches are connected, e.g. 1 single-pedal foot switch, 1 wireless double-pedal foot switch and 1 cabled double-pedal foot switch.

For detailed information on correctly connecting foot switches, refer to the instructions for use of the compatible foot switches.

For information on compatible foot switches refer to the chapter "15 Compatible equipment" on page 95.

4.6 Connection of peripheral equipment

Peripheral equipment can be connected to the LINK-IN socket and the LINK-OUT socket on the rear panel of the ESG-410. When connecting peripheral equipment, consider the following:

- The peripheral equipment must comply with respective IEC or ISO standards, e.g. IEC 60950-1/62368-1 for data processing equipment or IEC 60601-1 for medical equipment.
- For the connection the compatible communication cables MAJ-1871 (0.25 m) or MAJ-1872 (10 m) must be used.
- Also refer to the peripheral equipment's instruction for use.
- Both devices, the ESG-410 and the peripheral equipment, must be switched off before connecting.

4.7 Connection for automatic smoke evacuation


When the electrosurgical generator is combined with the Olympus high flow insufflation unit UHI-2, UHI-3 or UHI-4, smoke and mist can be evacuated simultaneously with the electrosurgical output.

4.7.1 Connecting the ESG-410 to UHI-2 or UHI-3

When connecting the ESG-410 to the UHI-2/3 also refer to the instructions for use of the UHI-2/3.

Preconditions

- UHI-2 or UHI-3 is available
- Communication cable MAJ-1872 (10 m) is available
- Adapter MAJ-1873 is available
- ESG-410 and UHI-2/3 are switched off

- 1) System connector
- 2) Adapter MAJ-1873 for UHI-2/3
- 3) Communication cable MAJ-1872

1. Ensure that both devices, the ESG-410 and the UHI-2/3 are switched off.
2. Connect the system plug of the adapter to the system connector of the UHI-2/3.
3. Fasten the screws on both sides of the adapter to secure the connection.


4. Firmly connect one plug of the communication cable to the LINK-IN socket or the LINK-OUT socket of the adapter.
5. Fasten the screws on both sides of the plug to secure the connection.
6. Firmly connect the other plug of the communication cable to the LINK-IN socket or the LINK-OUT socket of the ESG-410.
7. Fasten the screws on both sides of the plug to secure the connection.
8. Switch on both devices.

4.7.2 Connecting the ESG-410 to the UHI-4

When connecting the ESG-410 to the UHI-4 also refer to the instructions for use of the UHI-4.

Preconditions

- UHI-4 is available
- Communication cable MAJ-1872 (10 m) is available
- ESG-410 and UHI-4 are switched off

1. Ensure that both devices, the ESG-410 and the UHI-4 are switched off.
2. Firmly connect one plug of the communication cable to the LINK-IN socket or the LINK-OUT socket of the UHI-4.
3. Fasten the screws on both sides of the plug to secure the connection.
4. Firmly connect the other plug of the communication cable to the LINK-IN socket or the LINK-OUT socket of the ESG-410.
5. Fasten the screws on both sides of the plug to secure the connection.
6. Switch on both devices.

5 System settings

This chapter describes the correct operation of the touchscreen and the following basic settings are explained:

- Volume of activation tones
- Brightness of the touchscreen
- Alarm tones
- Touch tone
- Language selection
- Format for date and time
- Selection of monopolar modes
- Selection of startup settings

Furthermore, this chapter describes how to handle user-defined procedures and how to assign foot switches.

5.1 Operation of the touchscreen

Operating the touchscreen

- Only tap one button or pane at a time.
- Only tap the touchscreen from the front. Avoid oblique angles when tapping the touchscreen.

Do not use pointed objects, e.g. a pen, to tap the touchscreen.

Do not tap the touchscreen unintentionally.

Do not apply excessive loads to the surface of the touchscreen.

Remove dust and dirt by wiping the touchscreen with a cloth.

Operator's position

- The operator of the touchscreen must be in front of the electrosurgical generator.
- For parameter settings, the viewing distance must not exceed 50 cm.
- For watching the touchscreen, the viewing distance must not exceed 100 cm.
- All viewing angles, i.e. left, right, upper, lower angle, must not exceed 70°.

5.2 Settings – Setup

The following basic settings are selected on the [Setup] tab of the [Settings] screen.

- 1) Plus/minus buttons for volume
To increase/decrease the volume of the activation tone.
- 2) Plus/minus buttons for brightness
To increase/decrease the brightness of the touchscreen.
- 3) ON/OFF button for alarm tones
To switch the alarm tones for low and medium alarms on or off.
- 4) ON/OFF button for touch tone
To switch the touch tone on or off.
- 5) Options button for language
To select a preferred language.
- 6) Options button for date/time/format
To set the correct date and time.
- 7) Options button for monopolar modes
To select if the short list or the complete list of monopolar modes can be used.
- 8) Options button for startup settings
To select if the last HF output settings are restored after a startup of the ESG-410.
- 9) [Home] button
To return to the home screen.

5.3 Volume control of activation tones

The activation tones play an important role for the user to notice that the output of HF current is activated. If the activation tones are inaudible, then activation may not be recognized by the user.

- Adjust the volume level of the activation tones to the ambient noises. Make sure that the activation tones are at an audible volume level.

The volume of the activation tone can be adjusted in 2 ways:

- With the volume control on the back panel of the electrosurgical generator
- Via the touchscreen

During activation, it is impossible to enter the [Settings] screens. During activation, the volume must be adjusted by the volume control on the back panel.

The volume settings range from 1 to 10. The default setting is 7.

Setting the volume for the activation tones via the touchscreen

1. Tap [Settings] on the home screen.
 - ▶ The [Setup] screen is displayed.
2. In the volume pane tap [+] to increase the volume. Tap [–] to decrease the volume.

5.4 Brightness control

Use this function to control the brightness of the touchscreen. The brightness settings range from 1 to 10. The default setting is 5.

Setting the brightness

1. Tap [Settings] on the home screen.
 - ▶ The [Setup] screen is displayed.
2. In the brightness pane tap [+] to increase the brightness. Tap [–] to decrease the brightness.

5.5 Alarm tones

When an error message is displayed by the ESG-410 a certain alarm tone sounds, depending on the priority of the alarm. For alarms with a low or medium priority, the alarm tones can be switched off in 2 ways.

Switching off alarm tones on the [Settings] screen

1. Tap [Settings] on the home screen.
 - ▶ The [Setup] screen is displayed.
2. For [Alarm tone] tap the on/off button to turn the alarm tones on or off.
 - ▶ The selected setting is displayed [AUDIO ON] or [AUDIO OFF]

Switching off alarm tones on the error message

Every error message with a low or medium priority provides the function to switch the alarm tones on or off.

- Tap the button [AUDIO ON/OFF] to switch the alarm tones on or off.
 - ▶ This setting is saved for all low/medium priority alarms as long as the ESG-410 is switched on or as long as this setting is changed on any error message.

The default setting is [ON]. With every start of the ESG-410 the function Alarm Audio is automatically set to [ON].

Alarm tone for baby neutral electrodes

For safety reasons, the alarm tone for one low priority alarm cannot be switched off. This is the alarm tone when the current limit of 350 mA for baby neutral electrodes is reached. Refer to the section "6.3 Safety feature for baby neutral electrodes" on page 46.

5.6 Touch tone

Use this function to turn on or off the audible feedback when tapping the touchscreen. As a default setting, the audible feedback is turned off.

Switching the touch tone on/off

1. Tap [Settings] on the home screen.
 - ▶ The [Setup] screen is displayed.
2. Tap the on/off button to turn the audible feedback on or off.
 - ▶ The selected setting is displayed [ON or OFF]

5.7 Languages

Use the language function to select the preferred language. The texts on the touchscreen can be displayed in several different languages.

Selecting the language

1. Tap [Settings] on the home screen.
 - ▶ The [Setup] screen is displayed.
2. Tap the options button for [Language].
 - ▶ The list of available languages is displayed. The selected language is highlighted.
3. Tap the scroll arrows to scroll through the list.
4. Tap the preferred language.
 - ▶ A message box with a request for confirmation is displayed.
5. Tap [OK] to change the language setting.
 - ▶ The selected language is highlighted in the language list.
6. Tap [Back] to return to the [Settings] screen.
 - ▶ The selected language is displayed.

Scrolling through long lists

There are 3 ways to scroll through long lists:

- Tapping the scroll arrows.
- Moving the scroll box on the right hand side with the finger tip.
- Swiping the list with the finger tip.

5.8 Date, time and format

Use this function to set date and time of the system. The preferred date format is also set with this function.

Setting the format for date and time

1. Tap [Settings] on the home screen.
 - ▶ The [Setup] screen is displayed.
2. Tap the options button for [Date & time & format].
 - ▶ The [Date & time & format] screen is displayed.
3. Adjust the settings.
4. Tap [OK] to save the settings.
 - ▶ The display returns to [Settings] screen. The updated date/time/format settings are displayed.

5.9 Monopolar modes

The ESG-410 provides 2 sets of monopolar modes:

- A set with 4 standard monopolar modes
- A set with all 14 monopolar modes

Use this function to select the preferred set. The default setting is [All modes].

Selecting the set of monopolar modes

1. Tap [Settings] on the home screen.
 - ▶ The [Setup] screen is displayed.
2. Tap the options button for [Monopolar modes].
 - ▶ The [Monopolar modes] screen is displayed. The selected setting is highlighted.
3. Tap the button for the preferred setting.
4. Tap [Back] to return to the [Settings] screen.
 - ▶ The updated setting is displayed.

5.10 Startup settings

[Previous settings] and [All zero]

The ESG-410 provides 2 options on how to deal with HF output settings when the device is switched off:

- [Previous settings]
The selected HF output settings are stored. When the ESG-410 is switched on again, these settings are automatically displayed.
- [All zero]
The selected HF output and settings are deleted. When switching on the ESG-410 again, all HF output settings are automatically set to [Off].

The default setting is [All zero].

Selecting the settings for startup

1. Tap [Settings] on the home screen.
 - ▶ The [Setup] screen is displayed.
2. Tap the options button for [Startup settings].
 - ▶ The [Startup settings] screen is displayed. The selected setting is highlighted.
3. Tap the button for the preferred setting.
4. Tap [Back] to return to the [Settings] screen.
 - ▶ The updated setting is displayed.

[Last settings] button

If the ESG-410 is switched off unintentionally, then the last HF output settings can be restored by tapping this button.

If [All zero] was chosen as a startup setting, then this button is displayed on the top of the home screen when turning on the ESG-410.

The [Last settings] button is also available on the [Startup settings] screen.

5.11 Procedures

Use this function to save the current settings as a user-defined procedure or to select a user-defined procedure from a list.

Data that are saved within a user-defined procedure are:

- All selected modes
- All selected power settings
- All selected effect setting
- All selected level settings
- Instrument names
- The foot switch assignments
- The Autostart assignment
- The Autostart delay
- The setting for the baby neutral electrode
- The selection if an output socket is switched on or off.

User-defined procedures can be stored on a portable USB flash drive. Refer to the section “11.2.1 Data transfer” on page 87.

5.11.1 Saving a user-defined procedure

1. Select all settings on the touchscreen that are required for the new procedure.

2. Tap the [Procedure] button on the home screen.
 - ▶ The [Procedure] screen is displayed.
3. Tap the [Save] button.
 - ▶ The screen [New procedure] with an keyboard is displayed.
4. Enter a name for the new procedure.
5. Tap [OK].
 - ▶ The new user-defined procedure is saved.
 - ▶ The procedure name is displayed on the top right of the home screen.

Selecting a user defined-procedure

1. Tap the [Procedure] button on the home screen.
 - ▶ The [Procedure] screen is displayed.
2. Select a procedure name from the list.
 - ▶ The settings of the selected procedure are loaded.
 - ▶ The procedure name is displayed on the top right of the home screen.

Deviations from a user-defined procedure

If one or more parameters, e.g. the power level, mode name or effect, are manually changed while using a user-defined procedure, then an asterisk (*) is displayed in front of the currently selected procedure name.

Scrolling through long lists

There are 3 ways to scroll through long lists:

- Tapping the scroll arrows.
- Moving the scroll box on the right hand side with the finger tip.
- Swiping the list with the finger tip.

5.12 Assigning foot switches

There are 2 ways to assign foot switches to the output sockets:

- The initial assignment via the [Foot switch] button on the home screen.
- Changing the foot switch assignment clockwise with the black toggle button on the foot switch.

5.12.1 Assignment via the [Foot switch] button

Precondition

At least 1 foot switch is connected to the electrosurgical generator.

1. Tap the [Foot switch] button on the home screen.
 - ▶ The [Foot switch] screen is displayed
2. In the upper pane [Select socket], tap the desired socket.
 - ▶ The selected socket is highlighted.
3. In the lower pane [Assign foot switch], tap the desired foot switch.
 - ▶ The corresponding foot switch indicator is displayed within the selected socket.
 - ▶ The foot switch is assigned to the output socket.

5.12.2 Assignment via the toggle button on the foot switch

Precondition

At least 1 foot switch with a toggle button is assigned to one of the output sockets.

- Depress the black toggle button on the foot switch.
 - ▶ The foot switch indicator shifts clockwise to the next unoccupied output socket.
 - ▶ The foot switch is assigned to this output socket.

6 Safe use of neutral electrodes

Only for monopolar treatments, neutral electrodes are required to return the HF current from the active electrode of the HF instrument back to the electrosurgical generator. In bipolar electrosurgery both electrodes, the active and the return electrode, are implemented within the HF instrument. More explanations are described in the section "2.3.6 Monopolar electrosurgery" on page 15 of the safety chapter.

6.1 Split type and non-split type neutral electrodes and CQM

Contact quality monitoring (CQM) for split type neutral electrodes

If a split type neutral electrode is connected to the electrosurgical generator, CQM is able to detect 3 types of insufficiencies regarding the neutral electrode:

- The neutral electrode is detaching from the patient.
- The neutral electrode is disconnected, e.g. the plug is disconnected or the cable is torn from the neutral electrode.
- Short circuit of the neutral electrode, e.g. due to perspiration under the neutral electrode.

If one of these circumstances arises, then the following happens automatically:

- Activation of the electrosurgical generator stops and is prevented.
- A correspondent error message is displayed.
- An alarm tone sounds.
- The CQM indicator turns to red.

The correspondent error message instructs the user what to do to continue with the procedure.

No contact quality monitoring for non-split type neutral electrodes

When using non-split type neutral electrodes, then the contact quality monitor does not work. If the non-split type neutral electrode detaches from the patient, then no warning signal sounds. This can result in burns to the patient.

When connecting a non-split type neutral electrode, the following message box is displayed on the touchscreen:


The displayed message must be confirmed by the user for 2 reasons:

- To make the user aware that the risk of endogenous burns is high when using the non-split type of neutral electrodes.
- To avoid that a non-split type neutral electrode is used unintentionally.

6.2 Restriction to only split type neutral electrodes

As a safety feature, the ESG-410 provides a setting that prevents the use of non-split type neutral electrodes. If this setting is set to [Disable] and a non-split type neutral is (unintentionally) connected, then the following message box is displayed:

To switch this safety feature on/off proceed as follows:

1. Tap [Settings] on the home screen.
 - ▶ The [Settings] screen with the [Setup] tab is displayed.
2. Tap the [Maintenance] tab at the top.
 - ▶ The [Maintenance] tab is displayed.
3. Tap the [Non-split neutral electrode] button.
 - ▶ The [Neutral electrode] screen is displayed.
 - ▶ The current setting is highlighted.
4. Select the preferred setting.
 - ▶ The selected setting is highlighted and saved.
5. Tap [Home] to return to the home screen.

6.3 Safety feature for baby neutral electrodes

WARNING

Baby neutral electrode setting


When using neutral electrodes for newborns (patient weight less than 5 kg), the function [Baby neutral electrode] must be switched on. Otherwise, the electrosurgical generator can not detect the current flow's critical value of 350 mA. There is the risk of severe burns to the baby patient.

- Ensure that the function [Baby neutral electrode] is switched on, i.e. a baby symbol is shown within the neutral electrode indicator.
- Use the lowest power settings that are appropriate for the intended procedure.
- Use short activation intervals, as this reduces the heating between the neutral electrode and the newborn's skin to a minimum.

Functionality of the baby neutral electrode setting

A current flow of 350 mA serves as an indicator that the heating between the neutral electrode for newborns (patient weight less than 5 kg) and the baby's skin is getting to high. At this stage the user is warned by an alarm tone and the following error message:

If it is necessary to continue with the procedure, the user must observe the instructions that are stated within the error message:

- Minimize the activation time, e.g. by using short activation intervals and make pauses between the activation intervals.
- Reduce the power and effect settings or select another mode.


The error message erases when the user re-activates HF output or taps [OK]. However, the user still must pay special attention to the heating between the neutral electrode and the baby's skin.

As soon as the limit of 350 mA is reached again, the warning comes up again.

The default setting for the baby neutral electrode function is [OFF].

Setting of the baby neutral electrode

1. Tap the CQM info button on the bottom left of the home screen or on the set screen.
 - ▶ The [Neutral electrode setup] screen is displayed.
 - ▶ In the lower pane of this screen the [Baby neutral electrode] function is displayed.
2. Tap the options button to switch the function on or off.

If the function is set to [ON], then a baby icon is displayed within the CQM indicator. The selected setting remains, even if the electrosurgical generator is switched off and on again.

Data transfer of baby neutral electrode setting

The baby neutral electrode setting can be transferred to/from a portable USB flash drive only within a user-defined procedure. The baby neutral electrode setting is not stored on or loaded from the USB flash as a single setting, e.g. like the brightness or the touch tone.

When connecting a neutral electrode

If the baby neutral electrode function is set to [ON], then each time a neutral electrode is plugged to the ESG-410, the message box [Baby neutral electrode is selected] is displayed. This is to prevent that the warning message regarding excessive current comes up when working with neutral electrodes that are not for babies.

6.4 Conditions of the CQM indicator

The CQM indicator can adopt the following conditions:

No neutral electrode

There's no neutral electrode connected.

Indicator for split type neutral electrodes are green

- CQM detects sufficient contact between the neutral electrode and the patient's skin.
- Activation is possible.
- Continue with the procedure.

Indicator for split type neutral electrodes are red

- CQM detects an insufficient contact.

- Activation stops automatically and is prevented.
- An alarm tone sounds and a corresponding error message is displayed.
- Follow the instructions given in the error message.

Indicator for non-split type neutral electrodes is gray

- The connection of a non-split type neutral electrode is detected.
- Activation is possible.
- CQM does not work with non-split type neutral electrodes. There is a risk of burns to the patient if the neutral electrode detaches or gets damaged.

6.5 Selecting an appropriate neutral electrode

WARNING
Non-split type neutral electrodes

When using non-split type neutral electrodes, then the contact quality monitor does not work. If the non-split type neutral electrode detaches from the patient, then no warning signal sounds. This can result in burns to the patient.

Only use split type neutral electrodes. Refer to the section "15.2 Compatible neutral electrodes" on page 96.

Compatible neutral electrodes

The section "15.2 Compatible neutral electrodes" on page 96 lists split type neutral electrodes that are compatible for the use with this electrosurgical generator.

Only use compatible neutral electrodes as listed in the section "15.2 Compatible neutral electrodes" on page 96.

Only use self-adhesive neutral electrodes to prevent compression necrosis.

Do not use capacitive coupling neutral electrodes.

Immaculate neutral electrodes

- Check the use by date of the neutral electrode. Only use neutral electrodes with an unexpired use by date.
- Check the packaging of the neutral electrode. Only use neutral electrodes with an undamaged packaging.
- Check the neutral electrode for damages, modifications and sharp edges. Only use immaculate neutral electrodes.

Sufficient size

- Select the largest possible neutral electrode that is adequate for the patient. Consider the size of the neutral electrode particularly for obese patients. Consider the size of the neutral electrode particularly if high power levels are needed.

6.6 Connection of a neutral electrode to the electrosurgical generator


Specifications of the neutral electrode socket

The neutral electrode socket at the electrosurgical generator is suitable for plugs with the following specifications:

- Plugs with 2 receptacles
- 2.5 mm receptacle diameter
- Spacing: 10 mm


Baby neutral electrode function

If the baby neutral electrode function is set to [ON], then each time a neutral electrode is plugged to the ESG-410, the message box [Baby neutral electrode is selected] is displayed. This is to prevent that the warning message regarding excessive current comes up when working with neutral electrodes that are not for newborns.


6.6.1 Connecting a neutral electrode with a pre-attached cable

- Insert the plug of the neutral electrode into the neutral electrode socket on the front panel of the electrosurgical generator.

6.6.2 Connecting a neutral electrode without a pre-attached cable

To connect a neutral electrode without a pre-attached cable also refer to the instructions for use of the connecting cable and of the neutral electrode.

1. Lift the locking lever of the connecting cable.
2. Position the tab of the neutral electrode evenly between the clamp jaws of the connecting cable.
3. Completely push down the locking lever to lock the clamp.
4. Insert the plug of the connecting cable into the neutral electrode socket on the front panel of the electrosurgical generator.

6.7 CQM indicator verification

Verifying CQM when using a split type neutral electrode

- After attaching a neutral electrode to the patient, check the CQM indicator. The split type indicator must be green. Refer to the section "6.1 Split type and non-split type neutral electrodes and CQM" on page 44.

Verifying CQM when using a non-split type neutral electrode

- After attaching a neutral electrode to the patient, check the CQM indicator. The non-split type indicator must be gray. Refer to the section "6.1 Split type and non-split type neutral electrodes and CQM" on page 44.

If the indicator is red, green or gray although no neutral electrode is attached, there is a malfunction of the electrosurgical generator. In this case, do not use this electrosurgical generator.

7 Before use

This chapter describes all necessary steps to prepare the product for each use and provides detailed information on each mode.

7.1 Safety notes for before use

Before each operation

Any irregularity or damage of the electrosurgical generator and its equipment compromises the safety. Therefore, the inspections procedures as described in this chapter must be performed before every operation.

Do not use the electrosurgical generator if any irregularity or damage is suspected.

Do not use the equipment if any irregularity or damage is suspected.

- Prepare a spare electrosurgical generator or an alternative procedure to avoid interruption due to an unexpected electrosurgical generator failure during the procedure. Refer to the chapter "8.10 Troubleshooting" on page 78 to identify and correct failures.

Battery status of wireless foot switches

When the battery capacity of the wireless foot switch is insufficient, it is impossible to activate the HF output by pressing a foot switch pedal. When the foot switch unexpectedly fails during the procedure, there is the risk that a treatment, e.g. hemostasis, is interrupted.

Power levels and activation time

The selected power level and the activation time should be as low as possible to achieve the desired effect. If the power level initially is set too high or if the activation time is too long, then this can result in neuromuscular stimulation (violent spasms or muscle contractions) thermal injury to surrounding tissue, bleeding or inadequate treatment.

It needs to be considered that the minimum power level for the modes PulseCut Fast and PulseCut Slow is 20. The power levels of all other modes can be reduced to 1.

Apparent low output or failure of the HF surgical equipment

Apparent low output or failure of the HF surgical equipment to function correctly at the normal operating settings can indicate faulty application of the neutral electrode or poor contact in its connections. Therefore, before selecting a higher output power, the attachment of the neutral electrode and its connections must be checked.

HPCS not for conscious patients

Especially the modes that are supported by HPCS (High Power Cut Support) generate high current densities which can result in neuromuscular stimulation. This can cause pain and discomfort to patients without sedation, pain medication or general anesthesia.

Autostart for the mode BiSoftCoag

When using the Autostart function for the bipolar mode BiSoftCoag, then the HF output starts as soon as the HF instrument is in contact with tissue. Therefore, unintended tissue contact of the HF instrument can cause burns to the user and the patient.

The user should also be aware that the HF output is activated as long as the HF instrument is in contact with the tissue. The HF output only stops when the contact between the HF instrument and the tissue is released.

The power setting is limited to a maximum value of 50 when Autostart is switched on.

WARNING

Malfunction of pacemakers and defibrillators

Using HF equipment on patients with implanted electronic devices, e.g. cardiac pacemakers or cardioverter defibrillators, can cause failure of the implanted electronic device. Failure of the implanted electronic device will affect the heart and can result in cardiac arrest.

- Before the HF procedure, confirm its safety with a cardiologist or the manufacturer of the implanted electronic device.
- For monopolar procedures, position the neutral electrode so that the current pathway does not pass through or near the implanted electronic device and its lead system.
- Do not apply the HF instrument in close proximity to the implanted electronic device.

**WARNING****Flammable anaesthetics and gases**

If flammable gases, especially anaesthetics ignite, then they cause serious injuries to the patient, the user and the medical personnel.

- Take precautionary measures to keep away flammable gases from the site of intervention.
- Do not use flammable anaesthetics.

**WARNING****Flammable materials**

Sparks that occur during normal operation of the HF equipment can ignite flammable materials, e.g. absorbent cotton, gauze or surgical drapes and also body hair. This can result in serious injuries to the patient, the user and the medical personnel.

- Use soluble surgical lubricating jelly to cover hair close to the site of intervention.
- Do not apply materials like cotton or gauze simultaneously with the HF instrument to the site of intervention.

**WARNING****Irregularities and damage**

Any irregularity or damage of the electrosurgical generator and its equipment compromises the safety and can result in injury to the patient, the user and the medical personnel.

- Before every operation, inspect the electrosurgical generator as described in this chapter.
- Before every operation, inspect all equipment for irregularity or damage.
- Do not use the electrosurgical generator if any irregularity or damage is suspected.
- Do not use the equipment if any irregularity or damage is suspected.
- Refer to the section "8.10 Troubleshooting" on page 78.

**WARNING****Electrosurgical smoke**

Electrosurgical interventions produce smoke. When unfiltered smoke is inhaled, then this can affect the health of the user and the medical personnel.

- Wear high-filtration masks during the procedure.
- Make sure that a surgical smoke evacuator is used for the procedure.
- Make sure that adequate ventilation is provided.

**WARNING****Capacitive coupling due hybrid trocar systems**

When using monopolar active components then hybrid trocar systems, i.e. a combination of metal and plastic, can cause alternate site burns due to an accidental energy transfer. Only use all-metal or all-plastic trocar systems.

7.2 Inspection before use

7.2.1 Inspection regarding reprocessing

Make sure that the product has been properly reprocessed. Visually inspect the product thoroughly. The product must be visually clean.

Refer to the section "10 Reprocessing of the ESG-410" on page 84.

7.2.2 Testing the power supply

1. Switch on the electrosurgical generator.
2. Check that the power switch lights up.
3. Check that the touchscreen lights up.
 - ▶ The start screen with the Olympus logo is displayed shortly and the Olympus jingle is played.
4. Check that the home screen is displayed after the start screen.

7.2.3 Testing the touchscreen

1. Tap the pane of one of the output sockets, e.g. Monopolar 1.
 - ▶ The corresponding set screen is displayed.
2. Tap one mode button to change the mode setting.
3. Tap the plus and minus buttons.
4. Check that the power level setting changes accordingly.
5. Tap [Home] to return to the home screen.

7.2.4 Testing the foot switch

To properly inspect the foot switch, refer to the instructions for use of the foot switch that is used for the procedure.

7.2.5 Testing the alarm system

This section describes how to test the alarm system by exemplary testing the neutral electrode alarm.

1. Check that no neutral electrode is connected to the electrosurgical generator.
2. Check that the CQM indicator shows the status "No neutral electrode connected".
3. Try to activate any monopolar mode by pressing the corresponding foot switch pedal.
4. Check that the error message "E0202 Insufficient neutral electrode contact" is displayed and that an alarm tone sounds.
5. Check that the alarm tone stops and that the error message E0202 is closed after approximately 10 seconds.

7.3 Overview of modes

	Bipolar	Monopolar
Cutting	BipolarCut	Blend
	PlasmaCut	BlendCut
	PK LoopCut	FineCut
	PK MorceCut	PowerCut
	PK PureCut	PulseCut Fast
	PK SoftCut	PulseCut Slow
		Pure
Coagulation		PureCut
	AutoCoag	Fulgurate
	BiSoftCoag	ForcedCoag
	HardCoag	PowerCoag
	PK Coag	SoftCoag
	PK SoftCoag	Spray
	SalineCoag	SprayCoag

Additional modes within Olympus instruments

Special Olympus instruments, e.g. the POWERSEAL instruments, enable additional modes. When such an instrument is connected to the UNIVERSAL socket, the ESG-410 automatically recognizes the instrument, and the modes and settings are accessible. All required mode information of such instruments can be found in the corresponding instrument's instructions for use.

7.4 Description of modes

The following descriptions comply with the requirements of the standard IEC 60601-2-2.

¹⁾ The crest factor is specified at rated load.

²⁾ The power level is the continuous power (W) which can be applied during the activation of the electrosurgical generator. The actual applied power depends on the tissue characteristics, e.g. resistance.

7.4.1 Description of the bipolar cutting modes

BipolarCut

Fast cutting with small coagulation zone. High Power Cut Support (HPCS) enables fast start of cutting.

Control characteristics.....HPCS

Waveform Continuous sinusoidal, unmodulated

Crest factor¹⁾ 1.3 to 1.8

Range of power levels²⁾ 1 to 100 W

Effect levels 1 to 3

Maximum output voltage at effect level 1 425 V_p

Maximum output voltage at effect level 2 485 V_p

Maximum output voltage at effect level 3	645 V _p
Maximum output current	1.5 A
Rated load	500 Ω

PlasmaCut

Bipolar cutting in conductive fluid for transurethral resection (TURis).

Waveform	Continuous sinusoidal, unmodulated
Crest factor ¹⁾	1.3 to 1.8
Range of power levels ²⁾	1 to 320 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	597 V _p
Maximum output voltage at effect level 2	625 V _p
Maximum output voltage at effect level 3	700 V _p
Maximum output current	12 A
Rated load	75 Ω

PK LoopCut

Cutting of varying tissue structures, especially fibroid tissue with advanced bipolar technology. Intermittent output for fast and controlled cutting.

Waveform	Continuous sinusoidal, unmodulated
Crest factor ¹⁾	1.3 to 1.8
Maximum application time	60 s
Range of power levels ²⁾	1 to 200 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	440 V _p
Maximum output voltage at effect level 2	520 V _p
Maximum output voltage at effect level 3	590 V _p
Maximum output current	3 A
Rated load	200 Ω
Limited activation time	60 s

PK MorceCut

Cutting of varying tissue structures, especially fibroid tissue with advanced bipolar technology. Intermittent output for fast and controlled cutting and enhanced smoke reduction.

Waveform	Continuous sinusoidal, unmodulated
Crest factor ¹⁾	1.3 to 1.8
Maximum application time	60 s
Range of power levels ²⁾	1 to 300 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	500 V _p
Maximum output voltage at effect level 2	500 V _p
Maximum output voltage at effect level 3	560 V _p
Maximum output current	3 A
Rated load	75 Ω
Limited activation time	60 s

PK PureCut

Cutting of varying tissue structures using advanced bipolar technology. Fast Spark Monitor (FSM) ensures smooth and reproducible cutting in varying tissue, e.g. muscle and fat.

Waveform	Continuous sinusoidal, unmodulated
Crest factor ¹⁾	1.3 to 1.8
Maximum application time	60 s
Range of power levels ²⁾	1 to 200 W
Effect levels	1 to 3

Maximum output voltage at effect level 1	410 V _p
Maximum output voltage at effect level 2	440 V _p
Maximum output voltage at effect level 3	510 V _p
Maximum output current	1.5 A
Rated load	200 Ω

PK SoftCut

Cutting of varying tissue structures using advanced bipolar technology with enhanced smoke reduction. Intermittent output for fast and controlled cutting.

Waveform	Continuous sinusoidal, unmodulated
Crest factor ¹⁾	1.3 to 1.8
Maximum application time	60 s
Range of power levels ²⁾	1 to 200 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	430 V _p
Maximum output voltage at effect level 2	520 V _p
Maximum output voltage at effect level 3	610 V _p
Maximum output current	1.5 A
Rated load	200 Ω

7.4.2 Description of the bipolar coagulation modes**AutoCoag**

Slow bipolar coagulation for deep tissue effects. Automatic stop of activation avoids unintended thermal damage of surrounding tissue.

Waveform	Continuous sinusoidal, unmodulated
Crest factor ¹⁾	1.3 to 1.8
Range of power levels ²⁾	1 to 120 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	78 V _p
Maximum output voltage at effect level 2	134 V _p
Maximum output voltage at effect level 3	220 V _p
Maximum output current	1.4 A
Rated load	75 Ω

BiSoftCoag

Slow bipolar coagulation for deep tissue effects. Autostart option for activation when in contact with tissue.

Waveform	Continuous sinusoidal, unmodulated
Crest factor ¹⁾	1.3 to 1.8
Range of power levels ²⁾	1 to 120 W (max. 50 W with Autostart ON)
Effect levels	1 to 3
Maximum output voltage at effect level 1	95 V _p
Maximum output voltage at effect level 2	145 V _p
Maximum output voltage at effect level 3	220 V _p
Maximum output current	1.4 A
Rated load	75 Ω

HardCoag

Slow bipolar coagulation for deep tissue effects. Automatic stop of activation avoids unintended thermal damage of surrounding tissue.

Waveform	Continuous sinusoidal, unmodulated
Crest factor ¹⁾	1.3 to 1.8
Range of power levels ²⁾	1 to 120 W

Effect levels.....	1 to 3
Maximum output voltage at effect level 1	220 V _p
Maximum output voltage at effect level 2	220 V _p
Maximum output voltage at effect level 3	220 V _p
Maximum output current	3.3 A
Rated load	25 Ω

PK Coag

Coagulation of tissue with little sticking and carbonization using advanced bipolar technology and acoustic feedback on coagulation progress. Intermittent output for fast and controlled coagulation.

Waveform	Continuous sinusoidal, unmodulated
Crest factor ¹⁾	1.3 to 1.8
Range of power levels ²⁾	1 to 120 W
Effect levels.....	1 to 3
Maximum output voltage at effect level 1	150 V _p
Maximum output voltage at effect level 2	200 V _p
Maximum output voltage at effect level 3	220 V _p
Maximum output current	2.3 A
Rated load	75 Ω

PK SoftCoag

Coagulation of tissue with little sticking and carbonization using advanced bipolar technology with enhanced smoke reduction. Intermittent output for fast and controlled coagulation.

Waveform	Continuous sinusoidal, unmodulated
Crest factor ¹⁾	1.3 to 1.8
Range of power levels ²⁾	1 to 120 W
Effect levels.....	1 to 3
Maximum output voltage at effect level 1	150 V _p
Maximum output voltage at effect level 2	170 V _p
Maximum output voltage at effect level 3	230 V _p
Maximum output current	2.3 A
Rated load	75 Ω

SalineCoag

Bipolar coagulation in conductive fluid for transurethral resection (TURis).

Waveform	Continuous sinusoidal, unmodulated
Crest factor ¹⁾	1.3 to 1.8
Range of power levels ²⁾	1 to 200 W
Effect levels.....	1 to 3
Maximum output voltage at effect level 1	175 V _p
Maximum output voltage at effect level 2	230 V _p
Maximum output voltage at effect level 3	245 V _p
Maximum output current	5 A
Rated load	75 Ω

7.4.3 Description of the monopolar cutting modes**Blend**

Cutting with enlarged coagulation zone.

Waveform	Continuous sinusoidal, modulated
Crest factor ¹⁾	1.8 to 2.9
Range of power levels ²⁾	1 to 200 W
Effect levels.....	1

Maximum output voltage at effect level 1	2200 V _p
Maximum output current	1.0 A
Rated load	300 Ω

BlendCut

Cutting with enlarged coagulation zone.

Waveform	Continuous sinusoidal, modulated
Crest factor ¹⁾	3 to 4.5
Range of power levels ²⁾	1 to 200 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	935 V _p
Maximum output voltage at effect level 2	1140 V _p
Maximum output voltage at effect level 3	1400 V _p
Maximum output current	1.5 A
Rated load	500 Ω

FineCut

Fast cutting with small coagulation zone. High Power Cut Support (HPCS) enables fast start of cutting.

Control characteristics	HPCS
Waveform	Continuous sinusoidal, unmodulated
Crest factor ¹⁾	1.3 to 1.8
Range of power levels ²⁾	1 to 120 W
Effect levels	1
Maximum output voltage at effect level 1	560 V _p
Maximum output current	1.2 A
Rated load	500 Ω

PowerCut

Superficial coagulation with increased dissection capabilities.

Waveform	Continuous sinusoidal, modulated
Crest factor ¹⁾	3.5 to 6
Range of power levels ²⁾	1 to 120 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	1224 V _p
Maximum output voltage at effect level 2	1519 V _p
Maximum output voltage at effect level 3	2000 V _p
Maximum output current	1.2 A
Rated load	500 Ω

PulseCut Fast

Intermittent cutting and coagulation phases for controlled cutting. High Power Cut Support (HPCS) enables fast start of cutting. PulseCut Fast provides shorter coagulation phases than PulseCut Slow.

Control characteristics	HPCS
Waveform	Continuous sinusoidal, unmodulated
Crest factor ¹⁾	1.3 to 1.8
Range of power levels ²⁾	20 to 150 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	697 V _p
Maximum output voltage at effect level 2	706 V _p
Maximum output voltage at effect level 3	714 V _p
Maximum output current	1.5 A
Rated load	100 Ω

PulseCut Slow

Intermittent cutting and coagulation phases for controlled cutting. High Power Cut Support (HPCS) enables fast start of cutting. PulseCut Slow provides longer coagulation phases than PulseCut Fast.

Control characteristics.....	HPCS
Waveform	Continuous sinusoidal, unmodulated
Crest factor ¹⁾	1.3 to 1.8
Range of power levels ²⁾	20 to 150 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	671 V _p
Maximum output voltage at effect level 2	704 V _p
Maximum output voltage at effect level 3	714 V _p
Maximum output current	1.5 A
Rated load	100 Ω

Pure

Fast cutting with small coagulation zone.

Waveform	Continuous sinusoidal, unmodulated
Crest factor ¹⁾	1.3 to 1.8
Range of power levels ²⁾	1 to 300 W
Effect levels	1
Maximum output voltage at effect level 1	1300 V _p
Maximum output current	1.25 A
Rated load	300 Ω

PureCut

Fast cutting with small coagulation zone. High Power Cut Support (HPCS) enables fast start of cutting.

Control characteristics.....	HPCS
Waveform	Continuous sinusoidal, unmodulated
Crest factor ¹⁾	1.3 to 1.8
Range of power levels ²⁾	1 to 300 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	460 V _p
Maximum output voltage at effect level 2	525 V _p
Maximum output voltage at effect level 3	700 V _p
Maximum output current	1.5 A
Rated load	500 Ω

7.4.4 Description of the monopolar coagulation modes**Fulgurate**

Contact-free coagulation for fulguration and superficial tissue effect. Also used for cutting applications.

Waveform	Continuous sinusoidal, modulated
Crest factor ¹⁾	5 to 6.5
Range of power levels ²⁾	1 to 120 W
Effect levels	1
Maximum output voltage at effect level 1	3600 V _p
Maximum output current	1.5 A
Rated load	500 Ω

ForcedCoag

Coagulation for superficial tissue effect with little dissection capabilities.

Waveform	Continuous sinusoidal, modulated
Crest factor ¹⁾	5 to 6.5
Range of power levels ²⁾	1 to 120 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	1270 V _p
Maximum output voltage at effect level 2	1495 V _p
Maximum output voltage at effect level 3	2000 V _p
Maximum output current	1.2 A
Rated load	500 Ω

PowerCoag

Superficial coagulation with increased dissection capabilities.

Waveform	Continuous sinusoidal, modulated
Crest factor ¹⁾	3.5 to 6
Range of power levels ²⁾	1 to 120 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	1224 V _p
Maximum output voltage at effect level 2	1519 V _p
Maximum output voltage at effect level 3	2000 V _p
Maximum output current	1.2 A
Rated load	500 Ω

SoftCoag

Slow coagulation for deep tissue effects.

Waveform	Continuous sinusoidal, unmodulated
Crest factor ¹⁾	1.3 to 1.8
Range of power levels ²⁾	1 to 200 W
Effect levels	1 to 5
Maximum output voltage at effect level 1	150 V _p
Maximum output voltage at effect level 2	176 V _p
Maximum output voltage at effect level 3	198 V _p
Maximum output voltage at effect level 4	221 V _p
Maximum output voltage at effect level 5	235 V _p
Maximum output current	3.2 A
Rated load	50 Ω

Spray

Spark generation allows contact-free superficial coagulation.

Waveform	Continuous sinusoidal, modulated
Crest factor ¹⁾	7.5 to 9
Range of power levels ²⁾	1 to 120 W
Effect levels	1
Maximum output voltage at effect level 1	4300 V _p
Maximum output current	0.8 A
Rated load	500 Ω

SprayCoag

Spark generation allows contact-free superficial coagulation.

Waveform	Continuous sinusoidal, modulated
Crest factor ¹⁾	7.5 to 9
Range of power levels ²⁾	1 to 120 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	2520 V _p
Maximum output voltage at effect level 2	4300 V _p
Maximum output voltage at effect level 3	4300 V _p

Maximum output current 0.8 A
 Rated load 500 Ω

7.5 Tissue effects depending on the power level and effect level

7.5.1 Bipolar modes

Mode name	Results when increasing/decreasing the effect level	Results when increasing/decreasing the power level
BipolarCut PlasmaCut PK LoopCut PK MorceCut PK PureCut PK SoftCut	Increased / decreased thermal effect	Increased / decreased thermal effect and increased / decreased cutting capability
AutoCoag BiSoftCoag PK Coag PK SoftCoag	Increased / decreased coagulation speed	Increased / decreased coagulation depth for short time application. Nevertheless, to achieve the full potential coagulation depth a long application time with a low power setting is necessary.
HardCoag	Increased / decreased automatic end of procedure detection time	Increased / decreased thermal effect
SalineCoag	Increased / decreased thermal effect	Increased / decreased thermal effect

7.5.2 Monopolar modes

Mode name	Results when increasing/decreasing the effect level	Results when increasing/decreasing the power level
Blend BlendCut PulseCut Fast PulseCut Slow Pure PureCut	Increased / decreased thermal effect	Increased / decreased thermal effect and increased / decreased cutting capability
FineCut	n/a	Increased / decreased thermal effect and increased / decreased cutting capability
PowerCut	Increased / decreased dissection / coagulation	Increased / decreased thermal effect
Fulgurate PowerCoag	Increased / decreased coagulation / dissection capability	Increased / decreased thermal effect
ForcedCoag Spray SprayCoag	Increased / decreased coagulation speed	Increased / decreased thermal effect

Mode name	Results when increasing/decreasing the effect level	Results when increasing/decreasing the power level
SoftCoag	Increased / decreased coagulation speed	Increased / decreased coagulation depth for short time application. Nevertheless, to achieve the full potential coagulation depth a long application time with a low power setting is necessary.

8 Use

This chapter describes the following:

- Specific safety notes for use
- How to connect instruments
- How to set appropriate output settings for the intended procedure
- How to activate the output
- The alarm system with error messages and alarm tones
- Troubleshooting topics
- Overview of error messages

8.1 Safety notes for use

Maximum output voltage

When connecting instruments, avoid power settings (e.g. mode, effect), where the maximum output voltage exceeds the allowed peak voltage of the instrument. Refer to the section "7.4 Description of modes" on page 55. Also refer to the instructions for use of the instrument.

One instrument per output socket only

Do not connect more than 1 instrument per output socket.

No bipolar plugs with single connecting pins

If bipolar HF instruments are used that have single connecting pins without a fixed pin spacing, then there is the risk that the single pins are connected accidentally to the incorrect receptacles. This can result in damage of the electrosurgical generator.

Do not connect a 2-pin plugs with 4 mm pin diameters

Connecting a 2-pin plug with pin diameters of 4 mm (e.g. monopolar forceps with hand switch activation) to one of the MONOPOLAR output sockets can damage the electrosurgical generator.

WARNING

Patient and user injury due to hot tip

The HF instrument's tip remains hot enough for a few seconds to cause burns after the HF current is deactivated.

- Wait for a few seconds before laying the instrument aside.
- Do not place unused HF instruments on the patient.
- Store temporarily unused HF instruments in an electrically insulated container.

**WARNING****Connecting accessories**

Connecting accessories while the energy is on can result in injury and electric shock to the patient, the user and/or the medical personnel.
Ensure that the energy is off before connecting accessories.

**WARNING****Capacitive coupling with other surgical equipment**

An accidental energy transfer between the activated HF instrument and other surgical equipment can cause injury to the patient.
Only activate the HF instrument when it is in contact with target tissue.

**WARNING****Patient burns due to conductive fluids**

For monopolar treatments, conductive fluids like blood or saline that are in direct contact with or in close proximity to the active monopolar instrument can transmit electric current and heat to non-target tissue. This could cause unintended burns to the patient.
Before activating the instrument, aspirate any fluids from the treatment site.


**CAUTION****Power levels and activation time**

The selected power level and the activation time should be as low as possible for the intended purpose. If the power level initially is set too high or if the activation time is too long, then this can result in thermal injury to surrounding tissue, bleeding or inadequate treatment.

- For appropriate power levels, refer to the section “7.4 Description of modes” on page 55.

8.2 Description of the output sockets

8.2.1 BIPOLAR

The BIPOLAR socket is suitable for the following bipolar plugs:


- 2-pin plugs with 4 mm pin diameter, fixed pin spacing 28.6 mm
- 2-pin plugs with 4 mm pin diameter and a 3rd coding pin
- 2-pin plugs with 4 mm pin diameter, fixed pin spacing 22 mm, for foot switch activation.
- Coaxial plugs with 8/4 mm pin diameter, for foot switch activation.

Compatible Olympus instrument families

Bipolar HF cables (reusable or single-use) for:

- HICURA instruments
- HIQ+ instruments
- Everest instruments

8.2.2 MONOPOLAR 1 and MONOPOLAR 2


The MONOPOLAR sockets are suitable for HF instruments or HF cables with the following monopolar plugs:

- 3-pin plugs with 4 mm pin diameter
 - 1-pin plugs with 8 mm pin diameter
 - 1-pin plugs with 4 mm pin diameter
- Such a plug must only be connected to the right receptacle that is marked with a white ring.
- Coaxial plugs with 9/5 mm pin diameter

Compatible Olympus instrument families

- Monopolar electrosurgical pencils and instruments
- Monopolar HF cables (reusable or single-use) for:
 - HICURA instruments
 - HIQ+ instruments

8.2.3 UNIVERSAL 1 and UNIVERSAL 2

The UNIVERSAL sockets are suitable for special Olympus instruments or Olympus cables with instrument recognition.

Compatible Olympus instrument families

- POWERSEAL instruments
- PK instruments
- HF-resection electrodes supplied with single use HF cables

The reusable HF cable WA00014A is not compatible.

8.3 Connection of instruments

Before connecting

Only use this electrosurgical generator with compatible equipment. Refer to the chapter "15 Compatible equipment" on page 95.

- Check the maximum rated voltage of the other HF equipment. Make sure that the maximum output voltage of the electrosurgical generator does not exceed the rated voltage of the HF equipment that is used during the procedure.

Connecting

1. Refer to the instructions for use of the instrument.
2. Inspect the instrument, its connection cable and plug for scratches, cracks and other damages. Do not use the instrument if any kind of damage is observed.
3. Properly insert the plug of the instrument into the corresponding output socket on the front panel of the electrosurgical generator.

Limited usage time for HF resection electrodes with single use HF cables

The compatible HF-resection electrodes with single use HF cables have a limited usage time. 12 hours after the connected HF-resection electrode was activated for the first time, the electrosurgical generator automatically disables further activation. Corresponding messages are displayed on the touchscreen before the usage time expires.

8.4 Selection of appropriate output settings

Switching off output sockets

An output socket that is not needed can be switched off by tapping the toggle key on the top left of each output socket pane. If an output socket is switched off, unintended activation from this output socket is prevented. Furthermore, it improves the overview of the settings. Output socket panes that are switched off are darkened, so that only the relevant settings are displayed.

- 1) Toggle key to switch the display of the socket pane ON/OFF
- 2) Socket pane that is switched OFF
- 3) Socket pane that is switched ON

8.4.1 Mode selection

The mode is selected according to the type of procedure that will be performed and according to the instruments that will be used. Refer to the section "7.4 Description of modes" on page 55.

Selecting a mode

1. Tap the pane of the required output socket.
 - ▶ The set screen for the selected output socket is displayed.
2. Tap the mode name button.
 - ▶ A drop down list with modes names is displayed.
3. Tap the required mode name.
 - ▶ The mode is selected.

8.4.2 Power level

Selecting the power level

1. Tap the pane of the required output socket.
 - ▶ The set screen for the selected output socket is displayed.
2. Tap the plus or minus button to set the appropriate power level.

- For values up to 50, the power level increases/decreases in steps of 1.
- For values above 50, the power level increases/decreases in steps of 5.
- For values above 100, the power level increases/decreases in steps of 10.

To increase/decrease the power level fast, tap and hold the plus/minus button.

Actually applied power

The power levels that are displayed on the touchscreen represent the maximum possible output power in watts. The actually applied power depends on the tissue characteristics, e.g. resistance.

8.4.3 Effect level

The effect setting allows for customizing the mode settings and power settings. The user can influence the tissue effect by changing the effect level without increasing or decreasing the power setting.

Selecting the effect level

1. Tap the pane of the required output socket.
 - The set screen for the selected output socket is displayed.
2. Tap the arrow button to set the effect level.
 - The current setting is highlighted.
3. Tap the required effect level.
 - The select level is set.

8.5 Automatic instrument recognition


For UNIVERSAL sockets

Instruments that are connected to a UNIVERSAL socket are equipped with automatic instrument recognition. This function includes that appropriate output settings are set automatically. Depending on the connected instrument, the user is able to adapt these automatic settings within a restricted range.

Detailed information can be found in the corresponding instructions for use of the instrument.

8.5.1 The Seal mode for Olympus POWERSEAL instruments

This mode is available only when a corresponding Olympus POWERSEAL instrument is connected to the UNIVERSAL socket. The set screen for the Seal mode is configured differently compared to other modes.

- For the POWERSEAL instruments there are no cutting modes available. All cutting modes are set to [Off] automatically.
- The Seal mode is selected automatically.
- There are only 2 options for the power level: [1] or [–]
- Effect level settings are not available.

Messages regarding incomplete seal cycle

If a condition for an information message regarding an incomplete seal occurs, then the following takes place automatically:

- A 4-tone sequence sounds, indicating the incomplete seal.
- The output of HF current immediately stops.
- A corresponding information message is displayed, which provides information on remedial actions.

To clear the information message from the touchscreen, perform the remedial action that is displayed and Tap [OK] on the touchscreen. The information message can also be cleared by re-activating the power output with the hand switch or foot switch. For precise information on each information message regarding cause and remedy refer to the following table.

Info code	Message	Cause	Remedial action
i0763	Incomplete seal cycle 1. Inspect the seal. 2. Confirm no metal object is grasped. Remove excess fluids. 3. Reactivate instrument.	The electrosurgical generator detected a short circuit. The instrument jaws have contact with a metal part or there is excess fluid around the instrument jaws.	1. Release the hand switch or foot switch. 2. Open the instrument jaws and inspect the progress of the seal. 3. Remove excess fluids. 4. Confirm that no metal objects (staples, clips) are between the jaws. 5. Re-position the instrument to avoid metal contact. 6. Repeat the sealing process to complete the seal.
i0764	Incomplete seal cycle 1. Inspect the seal. 2. Confirm proper tissue contact. Grasp more tissue. 3. Reactivate instrument.	The instrument jaws are not closed completely. Insufficient amount of tissue is grasped.	1. Release the hand switch or foot switch. 2. Open the instrument jaws and inspect the progress of the seal. 3. If possible, re-grasp more tissue. 4. Close the instrument jaws completely. 5. Repeat the sealing process to complete the seal.
i0765	Incomplete seal cycle 1. Inspect the seal. 2. Hold button until complete seal cycle. 3. Reactivate instrument.	The seal cycle was interrupted before completion due to releasing the hand switch/ foot switch too early.	1. Release the hand switch or foot switch. 2. Open the instrument jaws and inspect the progress of the seal. 3. If the seal is incomplete, re-grasp the tissue (in the same position or overlapping the first seal a little). 4. Close the instrument jaws completely. 5. Repeat the sealing process to complete the seal.
i0766	Incomplete seal cycle 1. Inspect the seal. 2. Remove excess fluids. Clean and dry instrument jaws. 3. Reactivate instrument.	No tissue effect. The energy is absorbed by excess fluid or by coagulum that has built up at the instrument jaws.	1. Release the hand switch or foot switch. 2. Open the instrument jaws and inspect the progress of the seal. 3. Remove excess fluids. 4. Clean and dry the instrument jaws. 5. Re-grasp the tissue. 6. Repeat the sealing process to complete the seal.
i0767	Incomplete seal cycle 1. Inspect the seal. 2. Remove excess fluids. Clean and dry instrument jaws. 3. Reactivate instrument.	Automated generator time out. More time and energy is needed to complete the seal.	1. Release the hand switch or foot switch. 2. Open the instrument jaws and inspect the progress of the seal. 3. Remove excess fluids. 4. Clean and dry the instrument jaws. 5. Re-grasp the tissue. 6. Repeat the sealing process to complete the seal.

8.6 Selecting a user-defined procedure

For information on user-defined procedures refer to the section "5.11 Procedures" on page 41.

8.7 Activation of HF output

WARNING

Unintended tissue contact

When the output of the electrosurgical generator is active, then unintended contact between tissue and the active part of the HF instrument can cause burns to the patient, the user and the medical personnel.

- Store temporarily unused HF instruments in an electrically insulated container.
- Do not place unused HF instruments on the patient.

WARNING

Foot switch pedal or hand switch

Pressing the wrong foot pedal or hand switch can result in bleeding, perforation and burns to the patient.

- Make sure to press the correct foot switch or hand switch, i.e. yellow switch for cutting modes and blue switch for coagulation modes.
- Make sure that the required foot switch pedal is assigned to the correct output socket.

WARNING

Autostart

When using Autostart for the bipolar mode BiSoftCoag, then unintended tissue contact of the HF instrument can cause burns to the user and the patient.

- When using Autostart, be aware that the output of power starts as soon as the HF instrument is in contact with tissue.
- Avoid unintended tissue contact with the HF instrument.

WARNING

Visibility of the activated HF instrument

If the surgeon loses sight of the tip of the activated HF instrument during the endoscopic procedure, then this can result in bleeding, perforation and burns to the patient. Make sure that the tip of the activated HF instrument is always visible in the endoscopic image.

WARNING

Impeded performance

Coagulum or eschar that accumulates at the HF instrument can impede the performance of the electrosurgical generator. Increasing the power level can damage the electrosurgical generator and can result in injury to the patient.

- Carefully clean affected areas of the HF instrument. Refer to the instructions for use of the HF instrument.
- If necessary, replace the HF instrument.

CAUTION

Inappropriate settings

Inappropriate output of power can result in bleeding, perforation and burns to the patient. Before activating, confirm that the displayed settings are correct for the intended procedure.

**WARNING****Apparent low output**

Apparent low output or failure of the device to function correctly at the normal operating settings may indicate faulty application of the neutral electrode or poor contact in its connections.

Before increasing the intensity, check the adherence of the neutral electrode and its connections.

8.7.1 Considerations before activating**To prevent malfunction of electronic implants, e.g. pacemakers and defibrillators**

- Confirm with a cardiologist or with the manufacturer of the pacemaker/defibrillator that it is safe to perform the planned HF procedure.
- Do not apply the electrosurgical instrument in close proximity to the implanted electronic device.

To prevent procedural hazards and complications

- Before starting the procedure, check that all equipment is connected properly.
- Observe the instructions in the chapter "Inspection before operation".
- Prepare a spare electrosurgical generator or an alternative procedure to avoid interruption due to an unexpected electrosurgical generator failure during the procedure.

To prevent unintended current flows and HF leakage currents

- Ground the operating table.
- Place the patient on a dry, electrically insulating surface.
- The patient must not come in contact with earthed metal parts or with metal parts that have an appreciable capacitance to earth, e.g. the operating table.
- The patient's clothes must be dry.
- Avoid contact between the different skin surfaces (arms, legs) of the patient, e.g. place dry gauze between the body and arms and between the legs.
- Remove all metallic items from the patient, e.g. wristwatches, jewelry.
- Avoid skin contact between the patient and the user.
- Route all cables in a way that they are not in direct contact with the patient.
- Route all cables in a way that they are not in contact with other cables.
- Ensure that all equipment is connected properly.

8.8 Activating**WARNING****Baby neutral electrode setting**

When using neutral electrodes for newborn (patient weight less than 5 kg), the function [Baby neutral electrode] must be switched on. Otherwise, the electrosurgical generator can not detect the current flow's critical value of 350 mA. There is the risk of severe burns to the baby patient.

- Ensure that the function [Baby neutral electrode] is switched on, i.e. a baby symbol is shown within the neutral electrode indicator.
- Use the lowest power settings that are appropriate for the intended procedure.
- Use short activation intervals, as this reduces the heating between the neutral electrode and the baby's skin to a minimum.
- Refer to the section "6.3 Safety feature for baby neutral electrodes" on page 46.

Considerations before activating

- Confirm that the displayed settings are correct for the intended procedure.

- For monopolar treatments: Confirm that the CQM indicator for the neutral electrode is green (or grey if a non-split type neutral electrode is used).
- Make sure that the foot switch is assigned to the correct output socket.
- Make sure to use the correct foot switch pedal, i.e. the yellow pedal for cutting and the blue pedal for coagulation.
- Make sure to use the correct hand switch, i.e. the yellow cutting hand switch and the blue coagulation hand switch.

WARNING

Inadequate high output

Inadequate HF output, e.g. caused by a failure of the electrosurgical generator, can result in death or serious injury to the patient.

- Continuously watch the patient and the medical electrical equipment during use.
- If any abnormality of the patient or the medical electrical equipment is observed stop the application in a way that does not impair patient safety.
- Immediately stop the application of HF current as soon as inadequate HF output is observed.
- If the foot switch or the hand switch of the instrument does not react, then use the power switch on the ESG-410 as an emergency stop.

Considerations during activation

- Make sure that the tip of the activated instrument is always visible in the endoscopic image.
- Store temporarily unused instruments in an electrically insulated container. Do not place unused instruments on the patient.
- Carefully clean eschar affected areas of the instrument. Refer to the instructions for use of the instrument. If necessary, replace the instrument.

8.8.1 Functional principles of activating

Activating and stopping

Activation is possible only if the home screen or the set screen is displayed on the touchscreen.

To activate depress the corresponding pedal on the foot switch or press the corresponding hand switch of the instrument.

The HF output is activated as long as the foot switch pedal or the hand switch at the instrument is pressed. As soon as the corresponding switch is released, the HF output stops.

Activation tone


During activation an activation tone sounds.


It is impossible to change any setting during activation.

Background color

During activation the background color of the touchscreen adapts to the kind of output.


- For all cutting modes the background color changes to yellow.

- For all coagulation modes the background color changes to blue.

8.8.2 Automatic start for the bipolar mode BiSoftCoag

- 1) Autostart indicator
- 2) Options button to switch Autostart ON/OFF
- 3) Plus/Minus buttons for Autostart delay setting

If the tissue resistance is within a defined range, then the Autostart feature permits automatic activation of the output power as soon as both electrodes of the electrosurgical instrument touch the tissue. Consequently, activating by foot switch is not necessary. The HF output only stops when the contact between the HF instrument and the tissue is released.

It is possible to set a time delay so that the automatic activation does not start immediately but only after a certain time after the electrodes of the instrument get in contact with the tissue.

Autostart is available only for the bipolar mode BiSoftCoag.

8.8.3 Automatic stopping for the bipolar modes AutoCoag and HardCoag

When coagulating tissue the electrical resistance of the tissue increases. Tissue resistance is continuously measured by the electrosurgical generator so that it is able to detect when coagulation is completed. As soon as a limited value of tissue resistance is exceeded, the electrosurgical generator automatically stops the HF output.

The stopping of HF output is indicated by a special end tone.

Autostop is available only for the bipolar modes AutoCoag and HardCoag.

8.9 Alarm system

8.9.1 Alarm priorities

The electrosurgical generator is equipped with an intelligent alarm system. The alarm system determines alarm conditions on the base of multiple variables.

Alarm priorities

Depending on the potential risks, the alarms are classified in high, medium and low priority alarms. An alarm of higher priority overrides an existing alarm of lower priority. If more than one alarm condition of equal priority is determined, only the one that occurred first is reported.

The following table explains the 3 different conditions.


Alarm priority	Alarm tone	Caution symbol in error message	Meaning
High	3 short honks and 2 long honks for 2 cycles	The caution symbol is red and flashes.	- Serious system error. - The output is stopped automatically. - The electrosurgical generator automatically restarts. - If the restart runs successfully, then the user can continue to use this electrosurgical generator.
Medium	3 short honks	The caution symbol is yellow and flashes.	- With most medium priority alarms, the output stops automatically. - The user immediately must resolve the cause of the error as explained in the error message. - As soon as the error is resolved, the user can continue to use this ESG-410.
Low	2 short honks	The caution symbol is yellow without flashing.	- The user must take notice of the error message. - In some cases, the output is stopped automatically and the user must resolve the cause of the error as explained in the error message. - As soon as the error is resolved, the user can continue to use this electrosurgical generator.

Additional information on the alarm system

Additional information on the alarm system are described in the chapter "17 Additional information" on page 100.

8.9.2 Error messages

If an alarm condition occurs, then an alarm tone sounds and an error message is displayed on the touchscreen of the electrosurgical generator. The following figure shows an example:

- 1) Caution symbol and error code
- 2) Error title
- 3) Remedial actions
- 4) [OK] button

If the cause that released the alarm is resolved, then the error message is closed automatically after approximately 10 seconds. It is possible to tap [OK] in the error message to close the error message earlier.

If the cause of the alarm is not resolved, then the error message stays on the touchscreen. To clear the alarm, follow the instructions that are given in the error message.

8.10 Troubleshooting

8.10.1 General problems

Check the following table to identify and correct failures. If the problem cannot be resolved by the described remedial actions, stop using the electrosurgical generator and contact Olympus.

Be aware that repairs must only be performed by authorized service centers.

Problem	Cause	Remedy
The electrosurgical generator does not respond when pressing the power switch.	Improper connection of the power cable to the socket on the back panel of the electrosurgical generator or to the grounded mains electricity socket.	Check the power cable at the socket on the back panel and at the mains electricity socket for correct connection.
	The grounded fixed socket outlet has improper or no output voltage.	Check the grounded fixed socket outlet or use an alternative grounded fixed socket outlet.
	The power cable is damaged.	Check the power cable for damages and, if necessary, replace the power cable.
	Malfunction of the electrosurgical generator.	Do not use this electrosurgical generator. Contact Olympus.
After switching on the electrosurgical generator the touchscreen remains dark but sounds are audible.	Malfunction of the touchscreen.	Do not use this electrosurgical generator. Contact Olympus.
There is no output tone during activation.	The volume for the output tone is set to an inaudible level, e.g. due to high environmental noise.	Increase the volume either on the [Settings] screen or use the volume control on the back panel of the electrosurgical generator.
	Malfunction of the electrosurgical generator.	Do not use this electrosurgical generator. Contact Olympus.
The volume cannot be adjusted on the screen [Settings] or via the volume control on the back panel of the electrosurgical generator.	Be aware that only the volume for output tones is adjustable. Error-related alarm tones are not adjustable.	No action required.
	Malfunction of the electrosurgical generator.	Do not use this electrosurgical generator. Contact Olympus.
The electrosurgical generator does not respond to foot switch or hand switch activation.	Improper contact of the connection cable of the foot switch or the instrument to the corresponding socket of the electrosurgical generator.	Check the connection cable of the foot switch and the instrument for correct connection.
	The foot switch, the hand switch of the instrument or their connection cables are damaged.	Inspect the foot switch, the hand switch of the instrument and their connection cables for damages. If necessary, replace the foot switch, the instrument or the corresponding connection cable.
	The incorrect foot switch pedal or the incorrect hand switch is pressed.	Press the correct foot switch pedal or the correct hand switch of the instrument.
	The electrosurgical generator is not switched on.	Switch on the electrosurgical generator with the power switch.

Problem	Cause	Remedy
	2 switches for activating are pressed simultaneously, e.g. another foot switch pedal or hand switch of an instrument.	Release all switches for activation and check that no switch is pressed unintentionally. Then activate by using only 1 switch.
	A message box is displayed on the touchscreen.	Tap [OK] to close the message box or wait until the message box closes automatically after a few seconds.
	It is possible to activate only if the home screen or the set screen is displayed.	Return to the home screen or set screen.
	There is no mode selected, [Off] is displayed instead of a mode name.	Select an appropriate mode.
	The required output socket is switched off.	Tap the toggle key on the top right of the correspondent output socket pane to switch the output socket on.
	The power level is set to 0, [--] is displayed instead of a power setting.	On the set screen, tap the plus button to set an appropriate power setting.
	Malfunction of the electrosurgical generator.	Do not use this electrosurgical generator. Contact Olympus.
Foot switch pedal or hand switch of the instrument is pressed and the output tone sounds but no output is delivered.	Improper connection of the instrument plug with the output socket of the electrosurgical generator.	Check the plug of the instrument for correct connection.
	Malfunction of the electrosurgical generator.	Do not use this electrosurgical generator. Contact Olympus.
It is impossible to stop the output of the electrosurgical generator.	Malfunction of the foot switch or hand switch.	Immediately switch off the electrosurgical generator by pressing the power switch. Replace the foot switch or the instrument with hand switch.
	Malfunction of the electrosurgical generator.	Do not use this electrosurgical generator. Contact Olympus.
Its impossible to switch off the electrosurgical generator.	Malfunction of the electrosurgical generator.	Disconnect the power cable plug from the socket on the back panel of the electrosurgical generator and/or from the grounded mains electricity socket. Do not use this electrosurgical generator. Contact Olympus.

Problem	Cause	Remedy
The data transfer (Export) to the USB flash drive failed.	Other data are stored on the USB flash drive.	Ensure that the USB flash is completely empty.
	The USB flash drive is not connected properly or was removed too early.	Ensure that the USB flash drive remains connected until the export operation is completely finished. A corresponding message box is displayed on the touchscreen.
	The used USB flash drive is not compatible.	Compatible with USB 1.1. Use an USB flash drive with the following specifications: - max. 32 GB - FAT or FAT32 formatting (ExFAT, NTFS etc. are not supported)
The data transfer (Import) from the USB flash drive failed.	No data are stored on the USB flash drive for the import.	Check if the required data (user-defined procedures/setting) are available on the USB flash drive.
	The used USB flash drive is not compatible.	Compatible with USB 1.1. Use an USB flash drive with the following specifications: - max. 32 GB - FAT or FAT32 formatting (ExFAT, NTFS etc. are not supported)
The info message i1772 "Use of non-split neutral electrodes disabled" cannot be cleared by tapping the OK button.	A non-split neutral electrode is connected to the ESG-410 and the safety function regarding non-split neutral electrodes is set to [Disable].	1. Disconnect the neutral electrode from the ESG-410. 2. Change the setting for non-split neutral electrodes to [Enable]. Refer to the section "6.2 Restriction to only split type neutral electrodes" on page 45. Or: Connect a split-type neutral electrode so that the contact quality of neutral electrode is monitored during use.

8.10.2 Table of error messages

Most of the error messages that are displayed on the electrosurgical generator give a precise description about the cause of a problem and also how to resolve it. However, some of the error message are not self-explanatory. These error message are explained as follows:

Error code	Message	Cause	Remedial action
E####*	System error 1. The electrosurgical generator will automatically restart. 2. If the problem persists, contact Olympus.	Malfunction of the electrosurgical generator.	The electrosurgical generator automatically restarts. If the restart runs successfully, it is safe to use this electrosurgical generator. If any error occurs permanently or repetitively, do not use this electrosurgical generator and contact Olympus.
E2001	System error 1. Restart the electrosurgical generator. 2. If the problem persists, contact Olympus.	Malfunction of the electrosurgical generator.	Switch the electrosurgical generator off and then on again. If the restart runs successfully, it is safe to use this electrosurgical generator. If any error occurs permanently or repetitively, do not use this electrosurgical generator and contact Olympus.
E####*	System error 1. Press OK to continue. The electrosurgical generator is ready to use. 2. If the problem persists, restart the electrosurgical generator. 3. If the problem persists, contact Olympus.	Malfunction of the electrosurgical generator.	Tap the [OK] button to confirm the error message. If this does not solve the problem, switch the electrosurgical generator off and then on again. If the restart runs successfully, it is safe to use this electrosurgical generator. If any error occurs permanently or repetitively, do not use this electrosurgical generator and contact Olympus.

* "####" is the placeholder for various error codes.

9 After use

This chapter describes the necessary steps to disconnect the electrosurgical generator and all accessories.

Disconnecting

After use, disconnect the electrosurgical generator from the mains electricity and disconnect all accessories from the electrosurgical generator.

1. Switch off the electrosurgical generator by pressing the power switch.
 - The light of the power switch and the touchscreen turn off.
2. Disconnect all instruments from the corresponding output sockets.
3. Disconnect the connecting cable of the neutral electrode from the neutral electrode socket.
4. Disconnect the power cable from the mains electricity.
5. Disconnect the foot switches from the socket at the back panel of the electrosurgical generator.

10 Reprocessing of the ESG-410

After each use, immediately perform the procedures as described in the following sections. If cleaning is delayed, debris encrustation can become a source of infection. Encrustation can also result in malfunction of the electrosurgical generator.

Sockets

Cleaning the sockets of the electrosurgical generator can deform and corrode the contacts which will result in damage to the electrosurgical generator.

Do not clean the output sockets, the power socket or the foot switch sockets.

10.1 General information for reprocessing

The reprocessing cycle

Reprocess the product according to the instructions in this chapter.

Reprocess the product before each use.

Validated reprocessing methods

Validated means that the efficacy of the method or agent has been validated for reprocessing the product as described in this document.

Use of reprocessing agents

The manufacturer must have approved the reprocessing agent for cleaning and disinfection of medical devices and for the material to be reprocessed. Select contact times and concentrations of the reprocessing agent according to the instructions given by the manufacturer. The FDA, the EPA, the DGHM or a comparable institution must have approved and validated the reprocessing agent for its efficiency.

10.2 Cleaning

Suitable reprocessing agents for cleaning

Cleaning has been validated with an alkaline disinfectant based on low alcohol (<20%) and quaternary ammonium compounds (Sani-Cloth® Plus manufactured by PDI Professional).

Cleaning procedure

Perform the following cleaning procedure after use. If necessary, let the product cool down. If cleaning is delayed, debris will solidify and it can be difficult to effectively clean the surfaces.

1. Take a soft, lint-free cloth that is soaked with a suitable reprocessing agent.
2. Thoroughly wipe the surfaces of the product until visibly clean.
3. In case of heavy soiling, use additional cloths.
4. Before disinfection, let the surfaces dry or wipe them dry using a lint-free cloth.

10.3 Disinfection

Suitable disinfectants

Low-level disinfection has been validated with an alkaline disinfectant based on low alcohol (<20%) and quaternary ammonium compounds (Sani-Cloth® Plus manufactured by PDI Professional at 3 minutes contact time).

Inspection before disinfection

- Visually inspect the product thoroughly after cleaning. The product must be visually clean.
- If there are any signs of debris, repeat the cleaning process.

Disinfection procedure

1. Take a soft, lint-free cloth that has been soaked with a suitable disinfectant.
2. Thoroughly wet all surfaces of the electrosurgical generator by wiping with moderate pressure.
3. In a wet state, completely depress the power switch and wipe the power switch with moderate pressure.
4. In a wet state, completely rotate the volume control on the rear panel of the electrosurgical generator and wipe the volume control with moderate pressure.
5. The surfaces must remain visibly wet for the contact time that is designated by the disinfectant's manufacturer but not less than 3 minutes.
6. Use additional cloths to assure the designated continuous wet contact time.
7. Let the surfaces dry or wipe them dry using a lint-free cloth.

10.4 Reprocessing of other HF equipment

For the reprocessing or disposal of other HF equipment, e.g. neutral electrodes and HF instruments, refer to the corresponding instructions for use of these products.

11 Maintenance

This chapter provides information on the regular safety check and on the maintenance settings of the ESG-410.

Do not service or maintain any part of the product while in use with a patient.

Changing fuses

Be aware that fuses must only be changed by an authorized service center.

11.1 Regular safety checks

To ensure electrical safety and to maintain device performance, Olympus recommends to perform a safety check annually or according to national regulations and guidelines on all medical electrical equipment.

The regular safety check must only be performed by a qualified electrician with sufficient experience in maintaining medical electrical devices. Furthermore, the appropriate equipment for performing the safety check must be available. Details on the regular safety check are described in the maintenance manual for the ESG-410.

Alternatively, the electrosurgical generator and the foot switch can be sent to an authorized service center where the regular safety check is performed. Contact the local Olympus representative for information on authorized service centers. When sending the electrosurgical generator and the foot switch to an authorized service center, consider the information in the section "12.2 Shipment" on page 92.

11.2 Maintenance settings

The 5 buttons in the upper pane are addressed mainly to the physician and the medical personnel. Whereas the 7 buttons in the lower pane provide settings that are intended for the maintenance or service personnel.

The following subsections provide information on each function behind the buttons.

11.2.1 Data transfer

This function allows to export or import data to or from a USB flash drive that is connected to the rear of the ESG-410.

Export of settings

The currently selected settings, except for the Baby neutral electrode setting and the Autostart delay setting, are transferred to the USB flash drive. In detail these are:

- Language setting
- Volume level for activation tones
- Brightness level of touchscreen
- Touch tone setting
- Setting for date/time/format
- Set of monopolar modes
- Startup settings
- Non-split neutral electrode setting
- Maintenance reminder

Export of user-defined procedures

All user-defined procedures (including the baby neutral electrode setting and the Autostart delay setting) are transferred to the USB flash drive.

Baby neutral electrode and Autostart delay setting

For safety reasons, the setting [ON] for the baby neutral electrode function and also the Autostart delay setting are transferred to the USB flash drive only within a user-defined procedure. These 2 settings are not stored on or loaded from the USB flash drive as a single setting, e.g. like the brightness or the touch tone.

Import

- All settings are loaded from the USB flash drive.
- All user-defined procedures are loaded from the USB flash drive.

Be aware that all user-defined procedures that are currently stored on the ESG-410 are deleted, when importing settings from a USB flash drive.

Export error log

The internally stored error log data are transferred to the USB flash drive, e.g. for service purposes.

Requirements for the USB flash drive

The USB socket is compatible with USB 1.1.

- Storage capacity of max. 32 GB
- Formatting FAT32 or FAT (ExFAT, NTFS and others are not supported)
- Multiple partitions are not supported.
- Do not connect hard disks drives.

11.2.2 Error log

This function stores the log data of the last 1,000 errors. The last 100 errors are displayed.

Scrolling through long lists

There are 3 ways to scroll through long lists:

- Tapping the scroll arrows.
- Moving the scroll box on the right hand side with the finger tip.
- Swiping the list with the finger tip.

11.2.3 Event log

This function stores the last 10,000 events. The last 100 events are displayed.

11.2.4 Delete procedure

This function allows to delete user-defined procedures.

11.2.5 Non-split neutral electrode

This function allows to define if non-split type neutral electrodes can be used for monopolar procedures. For more information refer to the chapter “6 Safe use of neutral electrodes” on page 44.

11.2.6 Maintenance reminder

This function allows to set a date on which the electrosurgical generator reminds the user that the next regular safety check is due.

When the due date is reached then the maintenance reminder indicator is displayed on the top left of the touchscreen.

11.2.7 Safety test

This function is for purposes of the regular safety check only. This function is described in the maintenance manual of the electrosurgical generator.

11.2.8 System check

This function is for purposes of the regular safety check only. This function is described in the maintenance manual of the electrosurgical generator.

11.2.9 Software version and software update

Software

This function displays the software version that is currently installed on the electrosurgical generator.

Software update

The function [Software] also provides the button [Software update]. Software updates must only be performed by qualified personnel. To update the software, the new software version must be downloaded from the internet and stored on an USB flash drive. Follow the instructions that are displayed on the touchscreen.

Requirements for the USB flash drive

The USB socket is compatible with USB 1.1.

- Storage capacity of max. 32 GB
- Formatting FAT32 or FAT (ExFAT, NTFS and others are not supported)
- Multiple partitions are not supported.
- Do not connect hard disks drives.

11.2.10 Reset to defaults

This function resets all settings to factory defaults. All settings are deleted when using this function. User-defined procedures remain in the memory of the electrosurgical generator.

- Instrument settings: default
- Autostart delay: 1.0 s
- Non-split neutral electrode: Not allowed
- Monopolar modes: All modes
- Startup settings: All zero
- Maintenance Reminder: Disabled
- Language: English
- Touch tone: OFF
- Volume of activation tones: 7
- Brightness of the touchscreen: 5
- Date format: YYYY-MM-DD

11.2.11 PKCoag tones

For the mode PKCoag the 3 different activation tone sequences are available. This function allows to select one of the following:

- Ascending tone sequence.
- Descending tone sequence
- Constant tone sequence. This is the default setting.

For more information refer to the section "17.2.2 Activation tones for bipolar modes" on page 101.

This function is password protected. To receive the password contact Olympus.

11.2.12 Service

This function is for service purposes only.

12 Repair and shipment

Consider the following information if the ESG-410 needs to be repaired or shipped to Olympus.

12.1 Repair

Authorized service centers

Repairs may only be performed by qualified servicing personnel that has been authorized by Olympus. Contact an Olympus representative or an authorized service center for repair and warranty information.

Otherwise, Olympus can not be held responsible for the safety, reliability, and performance of the product.

Changing fuses

Be aware that fuses must only be changed by an authorized service center.

Loss of warranty

Any guarantee or warranty claims are forfeited if the user or an unauthorized servicing agency attempts repair of a product.

Description of defects

To enable the service center to perform repairs in a timely manner, send the product together with a detailed description of the damage or defect. The following particulars should be included:

- Catalog number
- Serial number or lot number (if possible)
- Precise description of the malfunction
- Delivery date
- Invoice copy (for possible guarantee or warranty claims)
- Internal order number of the customer (for correct accounting of the repair order)

12.2 Shipment

Contaminated products can cause infection

Sending used and contaminated products to Olympus represents a risk of infection. Olympus is entitled to refuse soiled or contaminated products for reasons of safety. Reprocess the product before sending in. If reprocessing will damage the product completely, clean the product as thoroughly as possible and mark it accordingly.

When sending products to an authorized service center for repair or warranty claims, use the original cardboard packaging for the transport of the product. If this is not possible, wrap each product individually in sufficient paper or sheets of foamed material and place them in a cardboard box. If applicable, use the respective transport protection (e.g. instrument tray) of the product.

Service centers do not accept warranty claims for damage caused by inadequate packaging.

13 Storage and disposal

This chapter provides information on the following:

- Storage of the electrosurgical generator
- Disposal of the electrosurgical generator

13.1 Storage

To prevent damage of the electrosurgical generator and the foot switch during storage, observe the following:

- Place the electrosurgical generator and the foot switch on a stable, level surface.
- Store the electrosurgical generator and the foot switch according to the environmental conditions described in the chapter "17 Technical data" on page 69.
- Do not store the electrosurgical generator and the foot switch in a location that is exposed to direct sunlight, x-rays, radioactivity or liquids.
- Do not store the electrosurgical generator and the foot switch in a location that is exposed to strong electromagnetic radiation, e.g. microwave medical treatment equipment, shortwave medical treatment equipment, magnetic resonance imaging equipment, radio or mobile phones.
- Do not subject the electrosurgical generator and the foot switch to strong impacts during storage.

13.2 Disposal

When disposing of the product or any of its components, follow all applicable national and local laws and guidelines. To prevent environmental contamination, reprocess the product before disposal.

Waste electrical and electronic equipment

In accordance with European Directive on waste electrical and electronic equipment (WEEE), the product must not be disposed of as unsorted municipal waste, but should be collected separately.

14 Ordering information

Description	Catalog number
Adapter MAJ-1873 for UHI-2/3	N3809530
Communication cable MAJ-1871, 0.25 m	N3809330
Communication cable MAJ-1872, 10 m	N3809430
Electrosurgical generator ESG-410	WA91307C
Energy cart TC-E410	K10030281
Foot switch, double pedal	WA91301W
Foot switch, double pedal, wireless	WA91302W
Foot switch, single pedal	WA91303W
HF cable MAJ-860, monopolar, 3.5 m, pin diameter: 8 mm	A02915A
HF cable MAJ-861, monopolar, 3.5 m, pin diameter: 4 mm	A02916A
High flow insufflation unit UHI-4	N3829660
Neutral electrode, for newborns, with cable, split type Patient weight < 5 kg	WA99512C*
Neutral electrode, for children, without cable, split type Patient weight 5 to 15 kg	WA99513A*
Neutral electrode, for children, with cable, split type Patient weight 5 to 15 kg	WA99513C*
Neutral electrode, for adults, without cable, split type Patient weight > 15 kg	WA99514A*
Neutral electrode, for adults, with cable, split type Patient weight > 15 kg	WA99514C*
Neutral electrode, universal size, without cable, split type Patient weight ≥ 5 kg (adults and children)	WA99516A*
Neutral electrode, universal size, with cable, split type Patient weight ≥ 5 kg (adults and children)	WA99516C*
Neutral electrode cable MAJ-814	E0427213
Neutral electrode cable	WA99510A*
Power cable type E/F	WA95621A
Power cable type B	WA95622A
Power cable type G	WA95623A

*Not licensed in accordance with Canadian law.


15 Compatible equipment

Olympus recommends the use of equipment as listed in this chapter. If combinations are used that are not listed in this chapter, the user takes the full responsibility. Future equipment can also be compatible. For more information, contact an Olympus representative.

Some of the products listed in this chapter may not be available in all sales territories.

15.1 System chart

Electromagnetic compatibility has been successfully tested with HF equipment that does not exceed a complete length (instrument including cables, adapters) of 4.5 m. Exceeding this length can cause electromagnetic disturbances.

- * Monopolar and bipolar HF instruments, also refer to section "8.2 Description of the output sockets" on page 65
- ** Neutral electrodes (with connecting cable), see below
- *** Neutral electrodes (without connecting cable), see below

- **** Electrosurgical instruments used with videoscopes and fiberscopes (without S-cord connector)
- 1) Not licensed in accordance with Canadian law.

15.2 Compatible neutral electrodes

15.2.1 Recommended compatible neutral electrodes

The following Olympus neutral electrodes are not licensed in accordance with Canadian law.

Olympus Catalog number	Description
WA99512C	Neutral electrode, for newborns, with cable Patient weight < 5 kg
WA99513A	Neutral electrode, for children, without cable Patient weight 5...15 kg
WA99513C	Neutral electrode, for children, with cable Patient weight 5...15 kg
WA99514A	Neutral electrode, for adults, without cable Patient weight > 15 kg
WA99514C	Neutral electrode, for adults, with cable Patient weight > 15 kg
WA99516A	Neutral electrode, universal size, without cable Patient weight ≥ 5 kg (adults and children)
WA99516C	Neutral electrode, universal size, with cable Patient weight ≥ 5 kg (adults and children)

15.2.2 Other compatible neutral electrodes

Manufacturer Model	Neutral electrode description
3M 8180F	Patient weight 5...15 kg
3M/Medline 9160, 9165	Patient weight 5...15 kg
Bowa 816-141*, 818-141*, 816-161*, 818-161*	Patient weight > 15 kg
Bowa 816-071*, 818-071*	Patient weight 5...15 kg
Bowa 816-072*, 818-072*, 816-092*, 818-092*, 816-112*, 818-112*	Patient weight > 5 kg

Manufacturer Model	Neutral electrode description
Conmed 51-7310, 51-7410	Patient weight > 15 kg
Conmed 425-2508, 425-2801, 51-7710	Patient weight 5...15 kg
Conmed 410-2000, 410-2100, 410-2200, 410-2400*	Patient weight > 2.2 kg
Covidien E7507, E7507-DB, E7509, E7509B	Patient weight > 13.6 kg
Covidien E7510-25, E7510-25DB*	Patient weight 2.7...13.6 kg

*Not licensed in accordance with Canadian law.

15.3 Compatible foot switches

Description	Catalog number
Foot switch, double pedal	WA91301W
Foot switch, double pedal, wireless	WA91302W
Foot switch, single pedal	WA91303W

16 Technical data

16.1 Specifications for the ESG-410

Size

Width	370 ± 2 mm
Height	241 ± 2 mm
Depth	501 ± 2 mm
Depth including connected cables	520 ± 10 mm
Weight	16.4 kg
Weight of packaging	2.48 kg

Power supply

Voltage range	AC 100 to 120 V/AC 220 to 240 V
Frequency	50 Hz/60 Hz
Maximum input power for AC 100 to 120 V	1150 VA
Maximum input power for AC 220 to 240 V	1500 VA
Fuse	2 x 250 V/T 10 AH

Classifications

According to IEC 60601-1, protection against electric shock:	
Medical electrical equipment	Class I
Applied parts	Type CF

Protection against harmful ingress of water	IP20
Mode of operation	Non-continuous operation

According to applicable European Union harmonization legislation	IIb
--	-----

Output

High frequency functions	Monopolar and bipolar
High frequency at rated load	344 to 516 kHz
Maximum HF power monopolar	300 W at 500 Ω
Maximum HF power bipolar	320 W at 75 Ω
All modes	25% duty cycle (e.g. activated for 10 s, deactivated for 30 s)

16.2 Ambient conditions for operation, storage and transport

Temperature

For operation and storage	+10 to +40 °C
For transport	-25 to +60 °C

Relative humidity

For operation and storage	30 to 85 %, non-condensing
For transport	10 to 85 %, non-condensing

Atmospheric pressure

For operation	700 to 1060 hPa (≤ 3000 m)
For storage and transport	500 to 1060 hPa (≤ 3000 m)

16.3 Plug requirements

For the MONOPOLAR output sockets

3-pin plugs.....	PinØ 4 mm
1-pin plugs.....	PinØ 8 mm
1-pin plugs.....	PinØ 4 mm
Coaxial plugs	PinØ 9/5 mm

For the BIPOLAR output socket

2-pin plugs.....	PinØ 4 mm, 28.6 mm pin spacing
2-pin plugs	PinØ 4 mm with 3rd coding pin
2-pin plugs	PinØ 4 mm, pin spacing 22 mm
Coaxial plugs.....	PinØ 4/8 mm

For the UNIVERSAL output sockets

Suitable for special Olympus instruments with instrument recognition

For the neutral electrode socket

2 pin plugs.....	Ø 2.5 mm, 10 mm pin spacing
------------------	-----------------------------

16.4 Resistance range for CQM

Allowable resistance range for the Contact Quality Monitor – CQM

For split type neutral electrodes	25 to 150 Ω
For non-split type neutral electrodes	<25 Ω

16.5 Specifications for cables

Power cables

WA95621A (for many European countries).....	Type E/F
WA95622A (for Canada, USA and other countries)	Type B
WA95623A (for United Kingdom and other countries)	Type G
Length	4.5 ±0.1 m
Plug on device side	angled

Communication cables

MAJ-1871	
Length	0.25 m
MAJ-1872	
Length	10 m

16.6 Specifications for the USB flash drive

Storage capacity	max. 32 GB
Formatting	FAT32 or FAT
ExFAT, NTFS and others are not supported.	
Compatible with USB 1.1	
Multiple partitions are not supported.	
Do not connect hard disks drives.	

17 Additional information

This chapter contains:

- Additional specifications for the alarm systems like log file and alarm priorities
- Information on alarm tones and activation tones
- Description of mode characteristics (acc. to IEC 60601-2-2)
- Output characteristic diagrams

17.1 Alarm system

An alarm condition is indicated about 1 second after detection and securing.

17.1.1 Log file for alarm conditions

The alarm system provides a log of occurrences of alarm conditions. The oldest log is overwritten by the latest log when the capacity of 1000 logs is reached.

If the alarm system is powered down, then the log is still maintained by an internal electrical power source. In case of a total loss of power, i.e. mains electricity and internal power source, the complete log is erased. The time that is needed for powering down is not captured within the error log.

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Regarding tone information for alarms, refer to the section below.

17.2 Tone information

17.2.1 Alarm tones

Depending on the priority of an error, the alarm tones sound different.

Low priority errors

Duty cycle..... Intermitting signal (number of burst = 2)
 Frequency 220/330 Hz
 Volume ≥ 65 dB(A)

Medium priority errors

Duty cycle..... Intermitting signal (number of burst = 3)
 Frequency 220/330 Hz
 Volume ≥ 65 dB(A)

High priority errors

Duty cycle..... Intermitting signal (number of burst = 10)
 Frequency 220/330 Hz
 Volume ≥ 65 dB(A)

- The volume of alarm tones is not adjustable.
- The alarm tones for errors with low or medium priority can be switched off. Refer to the section “5.5 Alarm tones” on page 38.

17.2.2 Activation tones for bipolar modes

Volume

With the lowest setting of the volume control, the volume of all activation tones is at least 40 dB(A).

Mode name	Tone sequence	Description
BipolarCut PlasmaCut PK PureCut PK SoftCut	Constant tone pitch at 554 Hz	Continuous tone during cutting.
PK LoopCut PK MorceCut	Constant tone pitch at 554 Hz	Continuous tone during cutting. In case of short circuit: Intermittent tone with 300 ms on time and 150 ms off time.
AutoCoag	Constant tone pitch at 440 Hz Alternating end tone at 880/1318 Hz	Continuous tone during coagulation. Alternating end tone when reaching a tissue resistance of $\geq 797 \Omega$. The activation automatically stops and the end tone sounds as long as the pedal/hand switch is pressed. If the initial resistance of the target tissue is high ($\geq 797 \Omega$) then the end tone immediately sounds. Activation is not possible. Use the mode BiSoftCoag instead.
BiSoftCoag PK SoftCoag SalineCoag	Constant tone pitch at 440 Hz	Continuous tone during coagulation.
HardCoag	Constant tone pitch at 440 Hz Alternating end tone at 880/1318 Hz	Continuous tone during coagulation. Alternating end tone when reaching a tissue resistance of $\geq 797 \Omega$. The activation automatically stops and the end tone sounds as long as the pedal/hand switch is pressed. If the initial resistance of the target tissue is high ($\geq 797 \Omega$) then the end tone sounds after the minimum application time of 3 – 5 seconds, depending on the selected effect level (1, 2 or 3).
PK Coag	3 different tone sequences to select from - Ascending - Descending - Constant	
PK Coag Ascending	5 ascending tone pitches, resistance dependent Amplitude modulated regrasp tone at 675 Hz	Continuous tone that changes its tone pitch stepwise in relation to the tissue resistance. The higher the tissue resistance, the higher the tone pitch. This tone sequence allows the user to monitor the coagulation progress if the visibility of the application area is limited by the instrument. Tone pitch 1: 440 Hz = tissue resistance from 11 to 214 Ω Tone pitch 2: 466 Hz = tissue resistance from 215 to 262 Ω Tone pitch 3: 493 Hz = tissue resistance from 327 to 406 Ω Tone pitch 4: 523 Hz = tissue resistance from 407 to 605 Ω Tone pitch 5: 554 Hz = tissue resistance above 605 Ω The regrasp tone sounds when the jaws of the HF instrument touch each other or if the tissue resistance is below 11 Ω . A corresponding message box is displayed on the touch screen.

Mode name	Tone sequence	Description
PK Coag Descending	5 descending tone pitches, resistance dependent Amplitude modulated regrasp tone at 675 Hz	Continuous tone with intermittent start tone. The tone sequence changes its tone pitch stepwise in relation to tissue resistance. The higher the tissue resistance, the lower the tone pitch. This tone sequence allows the user to monitor the coagulation progress if the visibility of the application area is limited by the instrument. Intermittent start tone: 351 Hz = tissue resistance from 11 to 274 Ohm Tone pitch 1: 227 Hz = tissue resistance from 275 to 346 Ohm Tone pitch 2: 173 Hz = tissue resistance from 445 to 551 Ohm Tone pitch 3: 135 Hz = tissue resistance from 552 to 799 Ohm Tone pitch 4: 122 Hz = tissue resistance above 800 Ohm The regrasp tone sounds when the jaws of the HF instrument touch each other or if the tissue resistance is below 11 Ω. A corresponding message box is displayed on the touch screen.
PK Coag Constant	Constant tone pitch at 440 Hz, with alternating tone pitches as end tone at 880/1318 Hz Amplitude modulated regrasp tone at 675 Hz	This is the default setting for the mode PK Coag. Continuous tone during coagulation. Alternating end tone when reaching a tissue resistance of $\geq 823 \Omega$. The activation does not stop automatically. As long as the pedal or hand switch is pressed, the residual activation is indicated by alternating tones. The regrasp tone sounds when the jaws of the HF instrument touch each other or if the tissue resistance is below 11 Ω. A corresponding message box is displayed on the touch screen.
Seal	Constant tone pitch at 440 Hz Intermittent end tone "complete seal cycle" at 880 Hz Alternating end tone "incomplete seal cycle" at 740/554 Hz	Continuous tone during vessel sealing. Intermittent end tone when the seal is complete. When the seal is incomplete, an alternating end tone sounds.

17.2.3 Activation tones for monopolar modes

Volume

With the lowest setting of the volume control, the volume of all activation tones is at least 40 dB(A).

Mode name	Tone sequence	Description
Blend BlendCut FineCut PowerCut Pure PureCut	Constant tone pitch at 554 Hz	Continuous tone during cutting.
PulseCut Fast PulseCut Slow	Alternating tone pitches at 554/440 Hz	Continuous tone during cutting.
ForcedCoag Fulgurate PowerCoag SoftCoag Spray SprayCoag	Constant tone pitch at 440 Hz	Continuous tone during coagulation.

17.2.4 Other tones

Startup tone (Olympus jingle)

Duration..... 0.375 s
 Frequency 500/1000/2000 Hz
 Volume ≥ 65 dB(A), not adjustable

Touch tone


Duration..... 0.09 s
 Frequency 500 Hz
 Volume ≥ 65 dB(A), not adjustable

17.3 Mode characteristics according to IEC 60601-2-2

Details on the mode characteristics are described in the section "7.4 Description of modes" on page 55.

17.4 Characteristics of High Power Cut Support (HPCS)

The following is an example of the characteristics of HPCS during PulseCut mode:




17.5 Output characteristic diagrams

The following diagrams show the HF output data for each mode.




Graph color:	Black	Red	Green	Blue	Purple
Effect setting:	1	2	3	4	5

17.5.1 Monopolar modes



Blend



BlendCut



FineCut



PowerCut




PulseCut Fast



PulseCut Slow




Pure



PureCut



Fulgurate



ForcedCoag



PowerCoag



SoftCoag

Spray



SprayCoag



17.5.2 Bipolar modes



BipolarCut



PlasmaCut



PK LoopCut



PK MorceCut



PK PureCut



PK SoftCut



AutoCoag



BiSoftCoag



HardCoag



PK Coag

PK SoftCosg

SalineCoag

18 Electromagnetic compatibility

Precautionary measures

Electromedical devices are subject to special precautionary measures concerning electromagnetic compatibility (EMC). This device is to be used only for the purposes described in this document and has to be installed, set up, and operated in compliance with the EMC guidelines.

WARNING

The use of accessories which are not approved by the manufacturer may result in an increase of electromagnetic emissions and the compliance with the stipulated limit values are not guaranteed anymore.

WARNING

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the electrosurgical generator, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

WARNING

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Guidelines and manufacturer's declaration – electromagnetic emissions		
This electrosurgical generator is intended for use in the electromagnetic environment specified below. The customer or the user of this electrosurgical generator should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	According to the intended purpose this electrosurgical generator emits electromagnetic energy. Electronic equipment close the electrosurgical generator may be affected, when HF output is activated. In stand-by mode HF energy is used for internal functions only. The RF emissions are very low. The requirements for group 1 devices have been tested and passed.
RF emissions CISPR 11	Class A	The electrosurgical generator is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is observed:
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	
		WARNING This electrosurgical generator is intended for use by healthcare professionals only. This electrosurgical generator may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating the device or shielding the location.

The essential performance requirement according to IEC 60601-1 is not applicable to the ESG-410. Basic safety is the performance requirement used during immunity testing.

Guidelines and manufacturer's declaration – electromagnetic immunity

This electrosurgical generator is intended for use in the electromagnetic environment specified below. The customer or the user of this electrosurgical generator should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U_T^* for ½ cycle (100% dip in U_T) 0% U_T for 1 cycle (100% dip in U_T) 0% U_T for 250 cycles (100% dip in U_T) 70% U_T for 25 cycles (30% dip in U_T)	0% U_T^* for ½ cycle (100% dip in U_T) 0% U_T for 1 cycle (100% dip in U_T) 0% U_T for 250 cycles (100% dip in U_T) 70% U_T for 25 cycles (30% dip in U_T)	Mains power quality should be that of a typical commercial or hospital environment. If the user of this electrosurgical generator requires continuous operation during power mains interruptions, it is recommended that the electrosurgical generator be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

*Note: U_T is the a.c. mains voltage before application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity			
This electrosurgical generator is intended for use in the electromagnetic environment specified below. The customer or the user of this electrosurgical generator should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V _{rms} and 6 V _{rms} in ISM bands ^c 150 kHz to 80 MHz	3 V _{rms} and 6 V _{rms} in ISM bands ^c	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b . Interference is possible in the proximity of devices bearing the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m	
	27 V/m 380 to 390 MHz	27 V/m	
	28 V/m 430 to 470 MHz	28 V/m	
	9 V/m 704 to 787 MHz	9 V/m	
	28 V/m 800 to 960 MHz	28 V/m	
	28 V/m 1.7 to 1.99 GHz	28 V/m	
	28 V/m 2.4 to 2.57 GHz	28 V/m	
	9 V/m 5.1 to 5.8 GHz	9 V/m	

Note 1: The higher frequency range applies for 80 MHz and 800 MHz.

Note 2: These guidelines may not be applicable in all situations. The spread of electromagnetic fields differs depending on the absorption and reflection of buildings, objects and people.

^a The field strength of stationary transmitters such as base stations of mobile phones and land mobile devices, ham radio stations, AM and FM radio and TV broadcasters are theoretically not 100% predictable. A site study is recommended to determine the electromagnetic environment regarding the stationary transmitters. If the field strength measured at the operation site of this device exceeds the compliance levels listed above, the device should be monitored to ensure that it functions as intended. Additional measures such as changing the orientation or operation site of the device may be required if unusual features of performance are observed.

^b The field strength should be less than 3 V/m for the frequency range of 150 kHz to 80 MHz.

^c The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz. The amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz and 50.0 MHz to 54.0 MHz.

19 Glossary

Ampere

The ampere or amp (symbol A) is the base unit of electric current in the International System of Units.

Applied part

An applied part, according to standard IEC 60601-1, is a part of a medical electrical equipment that comes into physical contact with the patient, e.g. an HF instrument. Or, it is that part of a medical electrical equipment through which current can flow between the patient and the medical electrical equipment, e.g. the neutral electrode socket.

Bipolar electrosurgery

Bipolar electrosurgical instruments are designed in a way that the current passes only between the 2 electrodes that are implemented within the same device, e.g. between the jaws of a bipolar forceps. In bipolar electrosurgery the current does not pass large parts of the patient's body but only through the tissue that is in contact with the two electrodes within the electrosurgical instrument. Only the target tissue is part of the electrical circuit. Compare with glossary entry "Monopolar electrosurgery".

Coagulation

A temperature of 60 to 70 °C in the area around the instrument's electrode gradually heats the tissue cells until their intracellular fluid escapes through the cell membranes. The results are protein denaturation, shrinking of the cells and of vessels, clotting of blood and hemostasis.

Coagulum

A semisolidified mass, as of blood or lymph, also called clot. Coagulum that sticks to the electrodes can impede the performance.

CQM

Abbreviation for Contact Quality Monitor. The CQM is able to detect detachment of the neutral electrode from the patient, if a split type neutral electrode is used. Refer to section "6.1 Split type and non-split type neutral electrodes and CQM" on page 44.

Current density

The greater the current is that passes through tissue, the greater the effect is on the tissue. With respect to electrode size, smaller electrodes provide a higher current density and result in a concentrated heating effect at the site of tissue contact. Following the same principle, the neutral electrode used in monopolar electrosurgery is large in relation to the electrode of the electrosurgical instrument. By dispersing the current that is returning to the electrosurgical generator the neutral electrode minimizes heat production at the attachment site.

Cutting

When the generated heat causes the intracellular fluid to rapidly evaporate, then this results in an explosive burst of the cell membrane. To achieve this, a very high voltage of approximately 200 Vp is required. Such a high voltage generates sparks between the tissue and the instrument's electrode/s. These sparks concentrate the HF current to specific points of the tissue that heats up to more than 100 °C. The fast repetition of this effect around the active electrode determines the cut.

Crest factor

The crest factor is a waveform parameter of alternating current. It shows the ratio of peak values to the effective value. High crest factors indicate peaks.

To give a reliable statement on the effect that the height of a crest factor has, it is necessary to distinguish if a cutting mode or a coagulation mode is used. For cutting modes it can be said, that the higher the crest factor is, the higher is the thermal effect on the tissue (including coagulation depth). For coagulation modes it can be said, that the higher the crest factor is, the stronger is the formation of sparks. This can be unwanted, as the risk of tissue carbonization rises with sparks. However, strong sparks can also be desired, e.g. for the SprayCoag mode, where very strong sparks are used for contact free, superficial coagulation.

DGHM

Abbreviation for "Deutsche Gesellschaft für Hygiene und Mikrobiologie", a German society for hygiene and microbiology which works to advance research in the fields of infectious diseases and microbiology.

Electrode

Electrodes are components in electrical engineering, which are used to make contact with a nonmetallic part within an electric circuit. The position of the electrodes determines the path that the current takes. Also refer to "Bipolar/Monopolar electrosurgery". In electrosurgery, the word "electrode" is often used as a synonym for an electrosurgical instrument, e.g. cutting electrode.

Electromagnetic compatibility

Electromagnetic compatibility (EMC) is the ability of electrical equipment and systems to function satisfactorily in their electromagnetic environment without introducing intolerable electromagnetic disturbances to anything in that environment, such as electromagnetic interference (EMI) or even physical damage in operational equipment.

Electrosurgery

Refer to glossary entries "Bipolar electrosurgery" and "Monopolar electrosurgery".

EMC

Abbreviation for electromagnetic compatibility. Refer to the glossary entry "electromagnetic compatibility".

EPA

Abbreviation for Environmental Protection Agency, an agency of the United States federal government.

Eschar

A dry crust produced by, e.g. a thermal burn, by a corrosive application or by gangrene. Eschar is high in resistance to electrical current. Keeping electrodes clean and free of eschar enhances the performance.

FDA

Abbreviation for Food and Drug Administration, a federal agency of the United States Department of Health and Human Services, one of the United States federal executive departments.

Grounding

Making an electric connection between a given point in a system and a local earth. The connection to local earth may be intentional or unintentional/accidental and it may be permanent or temporary. In electrical engineering "grounding" is colloquial used as a synonym for potential equalization.

HF

Abbreviation for high frequency. Refer to the glossary entry "High frequency".

HF instrument

Instrument that applies HF current for electrosurgical procedures, also called electrosurgical instrument.

HF output

Energy that is provided by the electrosurgical generator.

High frequency (HF)

Frequencies of alternating current above 100 kHz. High frequencies are used in electrosurgery, as lower frequencies provoke excessive neuromuscular stimulation and perhaps electrocution. High frequencies are also called radio frequencies (RF), as frequencies above 100 kHz can also transmit radio signals.

Leakage current

Leakage currents are current flows that pass another pathway than the initially intended electric circuit. A certain percentage of the current could pass via the patient and the operating table or the surgeon and then via the ground back to the electrosurgical generator. Earth leakage currents bear the risk of thermal tissue damages.

mA

Abbreviation for milliampere, refer to glossary entry "Ampere".

Monopolar electrosurgery

In monopolar electrosurgery the HF current passes through the patient as it completes the circuit from the active electrosurgical instrument to the neutral electrode and back to the electrosurgical generator.

Compare with glossary entry "Bipolar electrosurgery".

Radio frequency (RF)

Refer to the glossary entry "High frequency".

Rated load

The rated load is a value of non-reactive load resistance which, when connected, results in the maximum HF output power from each mode of the electrosurgical equipment.

Reprocessing

Reprocessing is the process that is carried out on a used medical device in order to allow its safe reuse including cleaning, disinfection and, if applicable, sterilization and related procedures. Reprocessing also includes the testing and restoring of the technical and functional safety of the used device.

Resistance

The electrical resistance is a fundamental property of a material or substance that quantifies how strongly it resists electrical current. The opposite of electrical resistance is the electrical conductance, that is the ease with which an electric current passes. Electrical resistance shares some conceptual parallels with the notion of mechanical friction. The standard unit of electrical resistance is the ohm (Ω), while electrical conductance is measured in siemens (S).

Spark discharge

A spark discharge is an electric discharge that is accompanied by a spark. It's a disruptive discharge.

Split type neutral electrode

With a split type neutral electrode, the conductive layer is split in 2 or more areas. If a part of such a neutral electrode detaches from the patient's skin, then the electric circuit is broken. This condition is detected by the contact quality monitor (CQM), which causes the electrosurgical generator to stop the output of HF current. For more information refer to the chapter "6 Safe use of neutral electrodes" on page 44.

UDI

Abbreviation for Unique Device Identifier. The Unique Device Identifier is a label on the medical device with human- and machine-readable information, which allows to uniquely identify the device through distribution and use.

Voltage

Voltage, also called electromotive force, is a quantitative expression of the potential difference in charge between two points in an electrical field. The greater the voltage is, the greater is the flow of electrical current through a conducting medium for a given resistance to the flow. The standard unit is the volt, symbolized by the uppercase letter V.

Waveform

In electronics, the term is usually applied to periodically varying voltages, currents or electromagnetic fields.

20 Index

A

Activation tones, for bipolar modes	101
Activation tones, for monopolar modes	103
Alarm priorities	77, 100
Alarm tones	101
All zero.....	41
Asterisk before procedure name	42
AutoCoag.....	57
AutoCoag, Autostop.....	76
Autostart	76
Autostart, BiSoftCoag.....	57
Autostop	76
Autostop, AutoCoag.....	57
Autostop, HardCoag	57

B

Baby neutral electrodes.....	46
Back panel.....	22
Basic settings	36
Bipolar coagulation modes, description.....	57
BipolarCut.....	55
Bipolar cutting modes, description	55
Bipolar electrosurgery	16
BiSoftCoag	57, 76
Blend	58
BlendCut.....	59

C

Capacitive coupling	17
Carbonization	15
Cardioverter defibrillator	16
CAUTION, signal word	8
Cleaning	84
Clinical benefits	10
Coagulation	14
Cochlear implants.....	16
Compatible equipment	95
Compatible foot switches.....	97
Compatible neutral electrodes.....	96
Conductive fluids	17
Contact quality monitoring (CQM)	44
Contraindications.....	9
Contraindications, relative	9
CQM – Contact Quality Monitoring.....	20
Cutting	14

D

Data transfer.....	87
--------------------	----

Delete procedure	88
Direct coupling	18
Disinfection	84

E

Earth leakage currents	17
Effects on tissue	62
Electric circuit	15
Electric current and biological tissue	13
Electromagnetic considerations for placement.....	30
Error log	88
Error message	25
Error messages	77
Error messages, table	81
Event log.....	88
Everest instruments	65
Export settings.....	87

F

FineCut	59
Fire and explosion	12
Flammable substances.....	12
Foot switch assignment	42
ForcedCoag.....	60
Front panel	21
FSM – Fast Spark Monitor.....	20
Fulgurate	60
Fuse changing	86
Fuses	92

H

HardCoag	57
HardCoag, Autostop	76
HF resection electrodes.....	66
HICURA instruments	65
HIQ+ instruments	65
Home screen	23
HPCS – High Power Cut Support.....	20

I

Icons on the touchscreen	28
Implants.....	16
Import settings.....	88
Incident reporting.....	10
Indications	9
Info message	25
Insulation failures.....	18
Intended use.....	9
Intended users.....	9

L

Last settings	41
Last values	41
Log file	100
Low output	52

M

Manufacturer	10
Message box	26
Metal implants	16
Mode selection	68
Mode tones	101
Monitoring equipment	17
Monopolar coagulation modes, description	60
Monopolar cutting modes, description	58
Monopolar electrosurgery	15
Monopolar modes setting	40

N

Needle monitoring electrodes	17
Neuromuscular stimulation	13, 52, 137
Neutral electrodes	16
Non-split neutral electrode setting	88
Non-split type neutral electrodes	44
NOTICE, signal word	8

O

Olympus jingle	104
----------------------	-----

P

Pacemakers	16
Patient target group	9
Physiological monitoring equipment	17
PK Coag	58
PKCoag tones	90
PK instruments	66
PK LoopCut	56
PK MorceCut	56
PK PureCut	56
PK SoftCoag	58
PK SoftCut	57
PlasmaCut	56
Power cables, compatible	94, 99
PowerCoag	61
PowerCut	59
POWERSEAL instruments	66
Procedures	41
PulseCut Fast	59
PulseCut Slow	60
Pure	60

PureCut	60
---------------	----

R

Rear panel	22
Regular safety checks	86
Reset to defaults	90
Resistances in tissue types	13

S

SalineCoag	58
Set screen	24
Smoke	12
SoftCoag	61
Software update	89
Software version	89
Split type neutral electrodes	44
Spray	61
SprayCoag	61
Startup settings	41
Startup tone	104
Symbols on the product/packaging	26

T

Toggle button of foot switch	42
Touch tone	104
Troubleshooting	78

U

UHI-2	32
UHI-3	32
UHI-4	34
Unintended current flows	17
USB flash drive, requirements	88, 89
User-defined procedure, deviations from	42
User-defined procedure, saving	41
User defined-procedure, selecting	42

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© Copyright 2021 Olympus Winter & Ibe GmbH

OLYMPUS DISTRIBUTION CENTERS

**Contact Information for
Customers**

www.olympus.com

United States of America and Canada

Olympus America Inc.

3500 Corporate Parkway, P.O. Box 610
Center Valley, PA 18034-0610
U.S.A.

Latin America

Olympus Latin America, Inc.

5301 Blue Lagoon Drive, Suite 290
Miami, FL 33126
U.S.A.

United Kingdom, Ireland

KeyMed Ltd.

KeyMed House, Stock Road
Southend-on-Sea, Essex SS2 5QH
United Kingdom

Europe

Olympus Europa SE & Co. KG

Amsinckstrasse 63
20097 Hamburg
Germany

Africa and Middle East

Olympus MEA FZ-LLC

Dubai Science Park, Nucleotide Complex,
Dubai, P.O Box 33607,
Registration No. 93456 (Dubai Creative
Clusters Authority),
United Arab Emirates

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Olympus Moscow Limited Liability Company

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OLYMPUS (BEIJING)

SALES & SERVICE CO., LTD.

A8F, Ping An International Finance Center,
No. 1-3, Xinyuan South Road, Chaoyang
District
Beijing, 100027, China

Olympus Hong Kong and China Limited

L43, Office Tower, Langham Place, 8 Argyle
Street,
Mongkok, Kowloon, Hong Kong

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Olympus Corporation

Shinjuku Monolith, 2-3-1 Nishishinjuku,
Shinjuku-ku, Tokyo 163-0914
Japan

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Asia, Australia and Oceania

Olympus Singapore Pte Ltd.

438B Alexandra Road #03-07/12
Alexandra Technopark Block B
Singapore 119968

Olympus AUSTRALIA PTY. LTD.

3 Acacia Place
Notting Hill, 3168
Melbourne
Building C Fern Tree Business Park

Olympus Medical Systems India Pvt Ltd

102-B, 1st Floor, Time Tower
M.G. Road, Gurgaon
Haryana-122001

Yuan Yu Industry Co., Ltd.

5 F., 37, Lane 365,
Yang-Guang St.,
Taipei 114, Taiwan

Published by

Olympus Surgical Technologies Europe

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
ESG-410

Sockets and screens Quick Reference Guide

ESG-410


Instrument sockets and compatible plug types

Sockets

Instrument plugs

MONOPOLAR


Monopolar instruments are not automatically recognized by the ESG-410.

⚠ WARNING


A 4 mm plug can only go into the white ring.

BIPOLAR


Bipolar instruments are not automatically recognized by the ESG-410.

UNIVERSAL


Only for Olympus instruments which are automatically recognized by the ESG-410.

Home screen

Set screen

Tap on a socket pane on the home screen (e.g. Bipolar) to adjust mode, power and effect for an instrument.

ESG-410

Foot switches, neutral electrodes and procedures

Quick Reference Guide


Foot switches

1. Connect a wired foot switch or a receiver to the foot switch sockets on the rear panel of the ESG-410.


Wireless

Wired

Rear panel of
ESG-410

Receiver

Available
foot switch types

Foot switch
indicators

Does the receiver match the foot switch?

There is one receiver for every foot switch.
The number stickers on both should be the same.

2. Tap the foot switch button on the home screen.
3. Tap the socket pane that requires a foot switch.
4. Once the selected socket is highlighted, tap the desired foot switch under [Select foot switch].

Foot switch

Select socket

Highlighted socket
pane

Unassigned
socket pane

Socket with an
assigned foot switch


Indicators of connected
foot switches

No foot switch
connected

Home

Neutral electrode setup screen

The [Neutral electrode setup] screen provides information on types and states of neutral electrodes.


Contact quality will be shown for split type/dual pad neutral electrodes. Contact quality is not known for non-split type/single pad neutral electrodes and won't be indicated.

Is the patient a baby?

Turn on baby mode to receive a warning message when the current flow through the neutral electrode pad is too high.

Procedure screen

To load a saved procedure, select a procedure from the list.

To save the current settings, tap [Save]. This will save:

- Modes, power, effect
- Foot switch assignment
- Autostart state
- Baby neutral electrode state

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ESG-410

Quick Reference Guide

	ESG-410	ESG-400	FT10	ForceTriad	ForceFX	System 5000	VIO 3	VIO 300 D	maXium	Mega Power	Autocon III 400	ARC 400
MONOPOLAR	Pure		Pure	Pure	Pure				PureCut	PureCut		
	Blend		Blend	Blend	Blend						C-Cut	
	PowerCut	PowerCoag		Valleylab		Pinpoint Coag	SwiftCoag	Swift Coag		Coag with Cut		Forced mixed
	PureCut	PureCut				Pure Cut	AutoCut	Auto Cut			PowerCut	Standard Laparoscopy
	BlendCut	BlendCut		Blend	Blend	Blend 1, 2, 3	DryCut	Dry Cut	Super Blend	Blend	C-Cut Lap-C-Cut	Dry
	PulseCut Slow	PulseCut Slow					EndoCut I/Q	Endo Cut I/Q	Polyp I, II		Papillo Cut Gastro-Cut	GastroLoop GastroKnife
	PulseCut Fast	PulseCut Fast					EndoCut I/Q	Endo Cut I/Q	Polyp I, II		Papillo Cut Gastro-Cut	GastroLoop GastroKnife
	FineCut	FineCut			Low			Precise Cut		ACE		Micro
	SoftCoag	SoftCoag	Soft				SoftCoag	Soft Coag			Standard Coag	Moderate
	ForcedCoag	ForcedCoag			Low (Dessicate)		ForcedCoag	Forced Coag	Forced Coag	Standard	Forced Coag	Forced non cutting
	SprayCoag	SprayCoag	Spray	Spray	High (Spray)	Spray Coag	SprayCoag	Spray Coag	Spray	Spray	Spray Coag	Spray
	Spray	SprayCoag	Spray	Spray	High (Spray)	Spray Coag	SprayCoag	Spray Coag	Spray	Spray	Spray Coag	Spray
BIPOLAR	PowerCoag	PowerCoag		Valleylab		Pinpoint Coag	SwiftCoag	Swift Coag		Coag with Cut		Forced mixed
	Fulgurate		Fulgurate	Fulgurate	Fulgurate	Standard Coag						
	BipolarCut	BipolarCut		Macro	Macro		AutoCut Bipolar	Bipolar Cut	Bipo Pure Cut		Bipolar Cut	Standard
	PlasmaCut	SalineCut	Bipolar Resection Cut				HighCut Bipolar	Bipolar Cut +	Uro Bip Pure Cut		Saline-C-Cut	Bipolar Resection (Cut)
	BiSoftCoag	BiSoftCoag		Low	Low (Precise)	Micro Coag	SoftCoag Bipolar	Bipolar Soft Coag	Macro Coag	Standard Micro Macro	Bip. Soft Coag Bip. Soft Coag with Autostop	Bipolar scissors
	AutoCoag	AutoCoag		Low & "A" ("A" = start & stop)			SoftCoag Bipolar & Auto Stop	Bipolar Soft Coag & Auto stop	Macro Auto Stop	Standard Micro Macro	Bivascular-Safe	Bipolar scissors
	SalineCoag	SalineCoag	Bipolar Resection Coag				SoftCoag Bipolar (Bipolar Resection)	Bipolar Soft Coag +	Bipol TUR Coag		Saline Coag	Bipolar Resection (Coag)
	HardCoag	HardCoag	LigaSure	LigaSure				BiClamp	SealSafe			Bipolar scissors

How to choose the right settings

3-step approach:

1. Find the equivalent mode in the mode lookup table.
2. Set the same power level that was used with the previous device.*
3. Generally, start with effect 2*.
(If small size electrodes are used, start with effect 1.)

ESG-410

ESG-400

FT10

ForceTriad

ForceFX

System 5000

VIO 3

VIO 300 D

maXium

Mega Power

Autocon III 400

ARC 400

* Except for PK modes or saline modes.

Start with the default settings that are automatically set, when the respective instrument/cable is connected.

Disclaimer:

The modes in this table are comparable on an application level (e.g. a pure cut mode without blend capability for smooth cutting with no or limited hemostasis). The mode lookup table should assist the user in selecting the ESG-410 mode which provides performance similar to modes experienced in a previously used

electrosurgical generator. Note that modes are never exactly the same between different generator manufacturers.

The mode lookup table is merely intended to be a guideline that the surgeon can assess for applicability. Depending on individual circumstances, surgeons may need to use other modes which are different from those provided in the lookup table.

OLYMPUS

Manufacturer:
Olympus Winter & Ibe GmbH
Kuehnstr. 61
22045 Hamburg, Germany

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		ESG-410	ESG-400	FT10	ForceTriad	ForceFX	System 5000	VIO 3	VIO 300 D	maXium	Mega Power	Autocon III 400	ARC 400
МОНОПОЛЯРНЫЕ	Резание	Pure		Pure	Pure	Pure				PureCut	PureCut		
		Blend		Blend	Blend	Blend						C-Cut	
		PowerCut	PowerCoag		Valleylab		Pinpoint Coag	SwiftCoag	Swift Coag		Coag with Cut		Forced mixed
		PureCut	PureCut				Pure Cut	AutoCut	Auto Cut			PowerCut	Standard Laparoscopy
		BlendCut	BlendCut		Blend	Blend	Blend 1, 2, 3	DryCut	Dry Cut	Super Blend	Blend	C-Cut Lap-C-Cut	Dry
		PulseCut Slow	PulseCut Slow					EndoCut I/Q	Endo Cut I/Q	Polyp I, II		Papillo Cut Gastro-Cut	GastroLoop GastroKnife
	Коагуляция	PulseCut Fast	PulseCut Fast					EndoCut I/Q	Endo Cut I/Q	Polyp I, II		Papillo Cut Gastro-Cut	GastroLoop GastroKnife
		FineCut	FineCut			Low			Precise Cut		ACE		Micro
		SoftCoag	SoftCoag	Soft				SoftCoag	Soft Coag			Standard Coag	Moderate
		ForcedCoag	ForcedCoag			Low (Dessicate)		ForcedCoag	Forced Coag	Forced Coag	Standard	Forced Coag	Forced non cutting
		SprayCoag	SprayCoag	Spray	Spray	High (Spray)	Spray Coag	SprayCoag	Spray Coag	Spray	Spray	Spray Coag	Spray
		Spray	SprayCoag	Spray	Spray	High (Spray)	Spray Coag	SprayCoag	Spray Coag	Spray	Spray	Spray Coag	Spray
БИПОЛЯРНЫЕ	Резание	PowerCoag	PowerCoag		Valleylab		Pinpoint Coag	SwiftCoag	Swift Coag		Coag with Cut		Forced mixed
		Fulgurate		Fulgurate	Fulgurate	Fulgurate	Standard Coag						
		BipolarCut	BipolarCut		Macro	Macro		AutoCut Bipolar	Bipolar Cut	Bipo Pure Cut		Bipolar Cut	Standard
		PlasmaCut	SalineCut	Bipolar Resection Cut				HighCut Bipolar	Bipolar Cut +	Uro Bip Pure Cut		Saline-C-Cut	Bipolar Resection (Cut)
	Коагуляция	BiSoftCoag	BiSoftCoag		Low	Low (Precise)	Micro Coag	SoftCoag Bipolar	Bipolar Soft Coag	Macro Coag	Standard Micro Macro	Bip. Soft Coag Bip. Soft Coag with Autostop	Bipolar scissors
		AutoCoag	AutoCoag		Low & "A" ("A" = start & stop)			SoftCoag Bipolar & Auto Stop	Bipolar Soft Coag & Auto stop	Macro Auto Stop	Standard Micro Macro	Bivascular-Safe	Bipolar scissors
		SalineCoag	SalineCoag	Bipolar Resection Coag				SoftCoag Bipolar (Bipolar Resection)	Bipolar Soft Coag +	Bipol TUR Coag		Saline Coag	Bipolar Resection (Coag)
		HardCoag	HardCoag	LigaSure	LigaSure				BiClamp	SealSafe			Bipolar scissors

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
WA91301W
WA91303W

Contents

1 General information	5
1.1 User instructions	5
1.2 Conventions throughout this document.....	5
1.3 Safety-related information in this document.....	5
1.3.1 General safety notes	6
1.3.2 Specific safety notes	6
1.3.3 Safety messages.....	6
2 Safety information.....	7
2.1 Intended use	7
2.2 User qualification	7
2.3 Environment of use	7
2.4 General safety notes	7
2.4.1 User-related error prevention	7
2.4.2 Risks regarding accessories	7
2.4.3 Risks regarding electric shock	8
2.4.4 Risks regarding fire and explosion	8
2.4.5 Risks regarding repair and maintenance	9
2.4.6 Risks regarding environmental conditions	9
3 General information	10
3.1 Maintenance reminder	10
3.2 Software version and software update	11
3.3 Reset to defaults	11
3.4 PK Coag tones	11
3.5 Overview of required tools	12
4 Inspection card.....	14
5 Visual inspection.....	20
5.1 Symbols and identification plate	20
5.2 Housing and fuse holder	21
5.3 UDI label	21
6 Electrical safety test	23
6.1 The safety test assistant	23
6.1.1 Starting the safety test assistant	24
6.2 Electrical safety	25
6.2.1 Checking the earth resistance and the earth leakage current	26
6.2.2 Checking the patient leakage current for HF output.....	26
6.3 Checking the DC resistance – Step 3	27
6.4 Measurements for the HF leakage current.....	28
6.5 Checking the HF leakage current for MONOPOLAR 1	29
6.5.1 Loaded condition.....	29
6.5.2 Unloaded condition	33
6.6 Checking the HF leakage current for MONOPOLAR 2	35
6.6.1 Loaded condition.....	35
6.6.2 Unloaded condition	37
6.7 Checking the HF leakage current for BIPOLAR.....	39
6.7.1 Loaded condition.....	39
6.8 Checking the HF leakage current for UNIVERSAL 1	41
6.8.1 Loaded condition.....	41

6.8.2 Unloaded condition	43
6.9 Checking the HF leakage current for UNIVERSAL 2	45
6.9.1 Unloaded condition - HF leakage current for UNIVERSAL 2	45
6.10 Checking the hand switch functions – Step 22	46
6.11 Checking the HF output power, open circuit and super ignition	47
6.11.1 Checking the HF output power for BIPOLAR – Steps 23-25	47
6.11.2 Checking the HF output power for MONOPOLAR 1 – Steps 26-30.....	48
6.11.3 Checking the HF output power for MONOPOLAR 2 – Steps 31-33.....	49
6.11.4 Checking the super ignition for UNIVERSAL 1 – Step 34	50
6.11.5 Checking the HF output power for UNIVERSAL 1 – Step 35.....	51
7 Contact quality monitoring (CQM).....	53
7.1 CQM test for bipolar modes	53
7.2 CQM test for monopolar modes	53
8 Final tests	54
8.1 Checking the power switch, sound volume and display	54
8.2 Enhanced system check	54
8.2.1 Testing without connections	54
8.2.2 Testing with connections	54
8.3 Foot switch function test.....	55
8.4 Communication test	55
8.5 USB socket test.....	56
8.6 Date and time	56
8.7 Pairing status of wireless foot switch	56
8.8 Inspection label	57
9 Disconnecting	58

1 General information

1.1 User instructions

This maintenance manual contains essential information on performing the regular safety check of the electrosurgical generator ESG-410 and its foot switch.

The complete set of instructions consists of the following documents:

- Instructions for use for the electrosurgical generator ESG-410
- Maintenance manual for the ESG-410 (this document)
- Instructions for use for all additional test equipment, e.g. the electrosurgical analyzer, the electrical safety tester.

Before use, thoroughly read this maintenance manual and the instructions for use of all equipment that will be used during the safety check procedures.

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Keep the maintenance manual in a safe, accessible location.

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	This is the safety alert symbol. It is used to alert the user to potential physical injury hazards. Observe all safety messages that follow this symbol to avoid possible injury.
	This symbol indicates additional helpful information.
1. 2. 3.	A numeration indicates a sequence of actions.
•	Bullet points indicate individual actions.
– – –	Dashes indicate the listing of data, options or objects.
1) 2) 3)	Numbers with right parenthesis name elements in illustrations.
[example]	Bracketed terms refer to elements in the graphical user interface or keys. Elements in the graphical user interface can be: - buttons - menu items - dialog elements

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DANGER

Indicates a hazardous situation which, if not avoided, will result in death or serious injury.

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2 Safety information

This chapter provides general safety information for regular safety check. This information is to be supplemented by the dangers, warnings and cautions given in the instructions for use of the electrosurgical generator, of the test equipment and in the instructions for use of any product being used with these products.

2.1 Intended use

Refer to the instructions for use of the electrosurgical generator.

2.2 User qualification

The regular safety check and software updates must only be performed by a qualified electrician with sufficient experience in maintaining medical electrical devices.

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2.4.1 User-related error prevention

WARNING
Improper use

The safety and effectiveness of electrosurgical interventions depend not only on the design of the equipment used, but also to a major extent on factors which are under the control of the user. It is therefore extremely important to read, understand and follow the instructions supplied with the electrosurgical generator and the accessories in order to ensure safety and effectiveness. Improper use will not only impede functions and prevent optimum performance, but can cause equipment damage and/or complications.

- Only use the electrosurgical generator as outlined in this maintenance manual.
- Before each use, inspect the equipment as outlined in this maintenance manual.

CAUTION
Regular safety checks

To maintain device performance and electrical safety, safety checks must be performed regularly on the electrosurgical generator and the foot switch. Otherwise, malfunction of the devices can occur which can result in injury to the patient.

- Olympus recommends to perform the safety check annually.
- Observe national statutory regulations.

2.4.2 Risks regarding accessories

CAUTION
Incompatible equipment and accessories

Using incompatible equipment and accessories can result in injury of the user as well as damage to the product.

- Only use compatible equipment and accessories. Refer to the chapter "Compatible equipment" in the respective instructions for use.
- Check the maximum rated voltage of all equipment and accessories.
- Make sure that the maximum output voltage of the electrosurgical generator does not exceed the lowest maximum rated voltage of any equipment and accessories that is used.

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The use of damaged equipment or of equipment with improper functioning can cause electric shock, mechanical injury, infection and thermal injury to the user. Do not use damaged equipment or equipment with improper functioning.

2.4.3 Risks regarding electric shock

WARNING

Grounding failure

There is a risk of electric shock if the housing of the electrosurgical generator is not grounded.

- Only connect the power cable to a properly grounded fixed socket outlet.
- Do not use an adapter for an incompatible grounded fixed socket outlet as it can impair safe operation.

WARNING

Precautions for electrical safety

Ignoring precautionary measures for electrical safety can result in electric shock to the user.

- Keep away liquids from all electrical equipment.
- Take appropriate precautions for electrical safety, e.g. use only isolated tools and equipment, follow the one hand rule.

CAUTION

Open housing

Opening the housing of the electrosurgical generator can result in electric shock to the user. Do not open the housing.

2.4.4 Risks regarding fire and explosion

WARNING

Flammable anaesthetics and fire supporting gases

The electrosurgical generator is not explosion proof. If flammable or fire supporting gases, especially anaesthetics ignite, then they cause serious injuries to the user.

- Take precautionary measures to keep away flammable and fire supporting gases from the site of intervention, e.g. nitrous oxide, oxygen.
- Do not use the electrosurgical generator in an explosive atmosphere.

WARNING

Flammable detergents and disinfectants

The igniting of flammable detergents and disinfectants can cause serious injuries to the user.

- If possible, only use non-flammable detergents and disinfectants.
- Before the electrosurgical generator is used, allow flammable substances to evaporate.

2.4.5 Risks regarding repair and maintenance

WARNING **Contamination**

Used medical products represent a risk of infection for the maintenance personnel. Before the maintenance work, clean and disinfect the products as stated in the respective instructions for use.

CAUTION **Repair**

Repairs must only be performed by authorized service centers. Do not attempt to modify or repair the electrosurgical generator or the foot switch. Contact the local Olympus representative or an distribution center for repair and warranty information.

2.4.6 Risks regarding environmental conditions

NOTICE

Environmental conditions

The use under improper environmental conditions can impede normal performance and/or can result in equipment damage.

Only use the complete equipment under the environmental conditions as described in the related instructions for use.

NOTICE

Interference with other equipment

The electrosurgical generator complies with the standards as described in the chapter "Electromagnetic compatibility" of the instructions for use. However, during activation, high frequency signals or spark discharge noise can disturb neighboring electrical equipment. Malfunction of the devices can occur.

- Follow the instructions regarding electromagnetic ambient conditions given in the chapter "Electromagnetic compatibility" of the instructions for use.
- Make sure that the devices are not used adjacent to or stacked with other equipment.
- Before use, thoroughly confirm the compatibility of all equipment.
- Do not use the electrosurgical generator in conjunction with:
 - Electrical equipment for which the safety against leakage current is not confirmed.
 - Electrosurgical equipment for which the safety in combined use is not confirmed.
- Do not loop cables.
- Do not bundle cables together with cables belonging to other medical equipment.
- Do not use the devices in a location that is exposed to strong electromagnetic radiation, e.g. microwave, short-wave medical treatment equipment, magnetic resonance imaging, radio, mobile phone.

3 General information

NOTICE

Inadvertent coupling and/or shunting of HF currents

- To avoid inadvertent coupling and/or shunting of HF currents around the resistor elements, keep the resistors at least 10 centimeters away from any metal surface including tabletops and other resistors. This is especially true if several resistors are connected in series or parallel to obtain a specified value.
- Do not allow the resistor bodies to touch each other.

Frequency of the safety check

To ensure electrical safety, Olympus recommends to perform the safety check annually or according to national regulations and guidelines on the following devices:

- The electrosurgical generator
- The foot switch

The safety check is mandatory after updating the software.

General requirements

- All steps must be performed in the order described.
- All tests must be done with fully functional and calibrated test equipment.
- Record the test results in the inspection card.
- Choose test cables with a length as short as possible, because cable inductance and stray capacitance can affect the readings.
- Carefully select suitable ground points to avoid ground loop errors in measurements.
- For tests that require a power cable, use the compatible power cable for the electrosurgical generator.
- Keep the signed inspection cards in a safe, accessible location.
- If any of the devices fail to meet any of the tests, contact Olympus.

Foot switches

Foot switches are tested together with the electrosurgical generator.

If a foot switch alone is tested, only the tests as described in the following sections are mandatory:

- "5.1 Symbols and identification plate" on page 20
- "8.2 Enhanced system check" on page 54 (with connections)
- "8.3 Foot switch function test" on page 55
- "8.7 Pairing status of wireless foot switch" on page 56 (if applicable for the foot switch)
- "8.8 Inspection label" on page 57 (if applicable for your country)

3.1 Maintenance reminder

This function allows to set a date on which the electrosurgical generator reminds the user that the next regular safety check is due. When the due date is reached then the maintenance reminder indicator is displayed on the top left of the touchscreen.

To set a due date for the next safety check proceed as follows:

1. Tap [Settings] on the home screen.
The [Settings] screen is displayed.
2. Tap [Maintenance].
The [Maintenance] tab is displayed.
3. Tap [Maintenance reminder] button.
The [Maintenance reminder] screen is displayed.
4. Tap the [Enable] button.

5. Select the settings for the due date of the next safety check.
6. Tap [Back] or [Home] to close the screen.
When the next safety check is due, the maintenance reminder icon will be displayed in the top left of the touchscreen.

3.2 Software version and software update

Observe the software version that is currently installed on the device.

1. Tap [Settings] on the home screen.
The [Settings] screen is displayed.
2. Tap [Maintenance].
The [Maintenance] tab is displayed.
3. Tap [Software Version].
The currently installed software version is displayed.

Software update

To update the software version tap the button [Software update] and follow the instructions that are displayed on the touchscreen.

3.3 Reset to defaults

This function resets all settings to factory defaults. User-defined settings are deleted when using this function. User-defined procedures remain in the memory of the electrosurgical generator.

- Instrument settings: default
- Autostart delay: 1.0 s
- Non-split neutral electrode: Not allowed
- Monopolar modes: All modes
- Startup settings: All zero
- Maintenance Reminder: Disabled
- Language: English
- Touch tone: Off
- Volume of activation tones: 7
- Brightness: 5
- Date format: YYYY-MM-DD

3.4 PK Coag tones

For the mode PKCoag 3 different activation tone sequences are available. This functions allows to select one of the following:

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For more information on the different tone sequences of these activation tones refer to the corresponding sections in the instructions for use.

To change the setting for the PKCoag activation tone proceed as follows:

1. Tap [Settings] on the home screen.
The [Settings] screen is displayed.
2. Tap [Maintenance].
The [Maintenance] tab is displayed.

3. Tap the [PK Coag tones] button.
A screen to enter a password is displayed.
4. Enter the password "Harmony".
The [PK Coag tones] screen is displayed.
5. Tap the arrow button and select a tone scheme according to your preferences.
6. Tap [Back] or [Home] to close the screen.

3.5 Overview of required tools

No.	Name	Model/REF	Specifications/Remarks
J001	Electrical safety tester	---	e.g.: Seculife ST, Unimet 1000 ST (Bender), QA-90 (Metron)
J002	Electrosurgical analyzer	---	e.g.: SECULIFE ES PRIME, Unitherm
J003	Digital multimeter	---	DC accuracy: < 1%, DC voltage range 500 V
J004	Load resistors (low power)	---	Low inductive part, short time load 5 Ω , 140 Ω , 170 Ω all 0.5 W (or more) and 1% tolerance or alternatively: resistor decade box
J005	Load resistor/s	WB979015	ESG Testbox or Low inductive part, 500 Ω (300 W, 5%)
J007	Cables, 4 mm to 4 mm (banana)	---	Connection cable: monopolar or bipolar output (1 x 4 mm) to various equipment (1 x 4 mm, banana plug), length 1 m 8 pieces required
J008	Cables, UNIVERSAL to 4 mm (banana)	WB979008	Connection cable: bipolar UNIVERSAL to various equipment (2 x 4 mm, banana plug), length 1 m 2 pieces required
J009	Cable, monopolar including hand switch	WB979013	Connection cable: monopolar output incl. hand switch to various equipment (1 x 4 mm, banana plug), length 3 to 4 m
J011	Cable, neutral electrode to 2 x 4 mm (banana)	WB979002	Connection cable: neutral electrode output (2 x 2.5 mm, 10 mm) to various equipment (2 x 4 mm, banana plug), length 1 m
J012	Communication cable	MAJ-1871	Length 0.25 m
J013	Crocodile clips (with 4 mm connection)	---	For connection to potential equalization bonding
J014	Foot switch	WA91301W WA91302W WA91303W	2 different foot switch types are required for the enhanced system check, i.e.: - 1 single-pedal foot switch (wired) and 1 double-pedal foot switch (wired) or - 1 single-pedal foot switch (wired) and 1 double-pedal foot switch (wireless) or - 1 double-pedal foot switch (wired) and 1 double-pedal foot switch (wireless)
J016	Power cable	---	Country specific

No.	Name	Model/REF	Specifications/Remarks
J021	USB flash drive	---	Storage capacity: max. 32 GB Formatting: FAT32 or FAT (ExFAT, NTFS and others are not supported) Multiple partitions are not supported. Do not connect hard disks drives
J032	2D barcode scanner	---	Capable to read GS1 data matrix e.g.: Keyence HR101 or a smartphone
J059	Hand switch adapter UNIVERSAL socket	WB979018	Connection cable: bipolar UNIVERSAL to various equipment (3 x 4 mm, banana plug)

4 Inspection card

- Copy this template.
- Record the values as described in the following tests.

Device description

Product name:	Electrosurgical generator: ESG-410	Foot switch/es:
REF number:		
Serial number (if available):		
Software version CPU no. and Package no.:		—
Maintenance Manual version (see back cover):		

Visual inspection

Test description	Requirement	Passed	
		Yes	No
Symbols, identification plate, product name	Present, legible and as outlined in this maintenance manual.		
Fuses	Match the indication on the identification plate.		
Housing	No damages		
2D barcode of UDI label	Barcode matches the human readable number.		
UDI label	Legible and properly attached.		

Earth resistance, earth leakage current and patient leakage current for HF

Safety test assistant step 2

Step no.	Test description	Requirement	Reading	Passed	
				Yes	No
2	Earth resistance	$\leq 0.3 \Omega$			
2	Earth leakage current Normal condition (NC)	$\leq 0.30 \text{ mA}$			
2	Earth leakage current Single fault condition (SFC)	$\leq 1.0 \text{ mA}$			
2	Patient leakage current for HF Normal condition (NC), AC	$\leq 10 \mu\text{A}$			
2	Patient leakage current for HF Normal condition (NC), DC	$\leq 10 \mu\text{A}$			
2	Patient leakage current for HF Single fault condition (SFC), AC	$\leq 50 \mu\text{A}$			
2	Patient leakage current for HF Single fault condition (SFC), DC	$\leq 50 \mu\text{A}$			

DC resistance

Safety test assistant step 3

Step no.	Test description	Requirement	Reading	Passed	
				Yes	No
3	DC resistance	$\geq 2 \text{ M}\Omega$			

Monopolar HF leakage current

Safety test assistant steps 4-13

The steps 6 to 9 are optional tests.

Step no.	Test description	Power	Effect	Load	Requirement	Measured value	Passed	
							Yes	No
4	Monopolar 1 Spray Coag (active electrode terminal)	120	3	500 Ω	$\leq 150 \text{ mA}$			
5	Monopolar 1 Spray Coag (neutral electrode terminal)	120	3	500 Ω	$\leq 150 \text{ mA}$			
6	Monopolar 1 ForcedCoag (active electrode terminal)	120	3	500 Ω	$\leq 150 \text{ mA}$			
7	Monopolar 1 ForcedCoag (neutral electrode terminal)	120	3	500 Ω	$\leq 150 \text{ mA}$			

Step no.	Test description	Power	Effect	Load	Requirement	Measured value	Passed	
							Yes	No
8	Monopolar 1 ForcedCoag (active electrode terminal)	120	3	open	≤ 150 mA			
9	Monopolar 1 ForcedCoag (neutral electrode terminal)	120	3	open	≤ 150 mA			
10	Monopolar 2 Fulgurate (active electrode terminal)	120	n/a	500 Ω	≤ 150 mA			
11	Monopolar 2 Fulgurate (neutral electrode terminal)	120	n/a	500 Ω	≤ 150 mA			
12	Monopolar 2 Fulgurate (active electrode terminal)	120	n/a	open	≤ 150 mA			
13	Monopolar 2 Fulgurate (neutral electrode terminal)	120	n/a	open	≤ 150 mA			

Bipolar HF leakage current

Safety test assistant steps 14-21

Step no.	Test description	Power	Effect	Load	Requirement	Measured value	Passed	
							Yes	No
14	BipolarCut (active electrode terminal 1)	100	3	500 Ω	≤ 70.7 mA			
15	BipolarCut (active electrode terminal 2)	100	3	500 Ω	≤ 70.7 mA			
16	Universal 1 PlasmaCut (active electrode terminal 1)	320	3	75 Ω	≤ 126.5 mA			
17	Universal 1 PlasmaCut (active electrode terminal 2)	320	3	75 Ω	≤ 126.5 mA			
18	Universal 1 PlasmaCut (active electrode terminal 1)	320	3	open	≤ 126.5 mA			
19	Universal 1 PlasmaCut (active electrode terminal 2)	320	3	open	≤ 126.5 mA			
20	Universal 2 PK PureCut (active electrode terminal 1)	200	3	open	≤ 100 mA			
21	Universal 2 PK PureCut (active electrode terminal 2)	200	3	open	≤ 100 mA			

Hand switch function

Safety test assistant step 22

Test no.	Test description	Requirement	Passed	
			Yes	No
22	Monopolar 1 and 2	Procedure passed		
22	Universal 1 and 2	Procedure passed		

HF output power, open circuit and super ignition tests

Safety test assistant steps 23-35

Step no.	Test description	Power	Effect	Load	Requirement	Measured value	Passed	
							Yes	No
23	HF output power Bipolar AutoCoag	120	3	75 Ω	$96.8 \text{ W} \leq P_{\text{out}} \leq 145.2 \text{ W}$			
24	HF output power Bipolar HardCoag	120	3	500 Ω	$7.44 \text{ W} \leq P_{\text{out}} \leq 11.1 \text{ W}$			
25	Open circuit test Bipolar HardCoag	120	1	5000 Ω	Message i0001 "Open circuit" is displayed			
26	HF output power Monopolar 1 SoftCoag	200	5	50 Ω	$158.4 \text{ W} \leq P_{\text{out}} \leq 237.6 \text{ W}$			
27	HF output power Monopolar 1 Spray Coag	120	3	500 Ω	$98.4 \text{ W} \leq P_{\text{out}} \leq 147.6 \text{ W}$			
28	HF output power Monopolar 1 PowerCoag	120	3	500 Ω	$96.0 \text{ W} \leq P_{\text{out}} \leq 144.0 \text{ W}$			
29	HF output power Monopolar 1 ForcedCoag	120	3	500 Ω	$94.4 \text{ W} \leq P_{\text{out}} \leq 141.6 \text{ W}$			
30	HF output power Monopolar 1 ForcedCoag	120	3	2000 Ω	$28.0 \text{ W} \leq P_{\text{out}} \leq 42.0 \text{ W}$			
31	HF output power Monopolar 2 Blend	200	n/a	300 Ω	$159.2 \text{ W} \leq P_{\text{out}} \leq 238.8 \text{ W}$			
32	HF output power Monopolar 2 Pure	300	n/a	300 Ω	$236.8 \text{ W} \leq P_{\text{out}} \leq 355.2 \text{ W}$			
33	HF output power Monopolar 2 Fulgurate	120	n/a	500 Ω	$95.2 \text{ W} \leq P_{\text{out}} \leq 142.8 \text{ W}$			

Step no.	Test description	Power	Effect	Load	Requirement	Measured value	Passed	
							Yes	No
34	Super Ignition Universal 1 Super ignition test	100	1	0 Ω	$2.4 \text{ A} \leq I_{\text{out}} \leq 3.5 \text{ A}$ and no error occurs			
35	HF output power Universal 1 PK PureCut	200	3	200 Ω	$153.6 \text{ W} \leq P_{\text{out}} \leq 230.4 \text{ W}$			

Contact quality monitoring (CQM)

Test description	Requirement	Passed	
		Yes	No
Any bipolar mode	CQM indicator "No neutral electrode connected" is displayed. Activation is possible.		
Any monopolar mode, standby, non-split type, 0 Ω , 5 Ω	Indicator for the non-split type neutral electrode is displayed in grey. Message box regarding non-split neutral electrode disabled is displayed.		
Any monopolar mode, standby, split type, 140 Ω	Indicator for the split type neutral electrode is displayed in green.		

Final tests

Test description	Requirement	Passed	
		Yes	No
Power switch LED, sound volume and display	LED shines green Volume and brightness adjustable		
Enhanced system check (with and without connections)	Procedure passed		
Foot switch function test	Procedure passed		
Communication test (LINK-IN/LINK-OUT sockets)	Error message [E0394 Communication error] is displayed.		
USB socket	Data transfer successfully		
Date and time	Correct date and time displayed		
If applicable: Pairing status of wireless foot switch	Pairing labels match. Connection indicator of receiver shines green when activating.		
If applicable: Inspection label	New inspection label attached with due date.		
Due date for next safety check:			

Result and signature

The device/s as stated under "Device description" meet/s the test criteria:	Yes	No
Date:		
Name of testing person:		
Test organization:		
Signature:		

5 Visual inspection


Preconditions

- The power switch of the electrosurgical generator is not depressed.
- The electrosurgical generator is not connected to the mains electricity.

5.1 Symbols and identification plate

Electrosurgical generator ESG-410

1. Check that all symbols and markings on the front and rear panel are present and legible. Refer to the instructions for use.
2. Check that the product name on the front panel is clearly visible.


3. Check that the identification plate of the ESG-410 shows details as outlined in the figure above, including the corresponding REF number.
4. Record the test results.

Foot switches

1. Check that the foot switch does not have:
 - any dents, cracks, kinks, or deformations
 - any deep scratches
 - any corrosion
 - any missing or loose parts
2. Check that all labels are present and legible.
3. Check that the identification plate of the foot switch is present and legible.
4. Record the test results.

5.2 Housing and fuse holder


1) Power entry fuse holder

1. Check the housing, the front panel and the rear panel for damages.
2. Check that the fuses match the indication on the identification plate. The arrows in the figure above show the position of the fuse labels.
3. Record the test results.

5.3 UDI label

The following figure exemplarily shows the following:

- Top: exemplary display of the scanning tool
- Bottom: exemplary corresponding UDI label on the device

- 1) Position number
- 2) GTIN-14
- 3) Position number
- 4) Production date
- 5) Position number
- 6) Serial number

Original GTINs for the ESG-410

GTIN for ESG-410 WA91307C: 04042761086349
 GTIN for ESG-410 WA91307J: 04042761086325
 GTIN for ESG-410 WA91307W: 04042761086332
 GTIN for ESG-410 WA91317J: (to be determined)

1. Check that the UDI label is legible.

2. Check that the UDI label is properly attached in the recess on the front panel.
3. Scan the 2D barcode with the tool J032.
4. Check that displayed numbers are equal to the human readable numbers of the UDI label.
5. Check that displayed numbers are equal to the original GTIN.
6. Record the test results.

6 Electrical safety test

NOTICE


If a 4 mm plug is incorrectly connected to the MONOPOLAR socket, this socket can be destroyed during activation.

Only connect a 4 mm plug to the receptacle on the right side.

6.1 The safety test assistant

The safety test assistant provides an easy approach to the required electrical measurements.

- 1) Step number 20/35
- 2) Progress bar with 35 single step indicators
- 3) Pane with test description
- 4) Pane with parameters
Displays the currents settings for this test and, if applicable, the required load resistance.
- 5) [Step done] button
To finalize the current step of the electrical safety test.

- 6) [Next] button
To go to the next step of the electrical safety test.
- 7) [Home] button
To go to the home screen.
- 8) [Back] button
To go back to [Settings] screen.
- 9) [Previous] button
To go to the previous step of the electrical safety test.

Progress bar and records

The progress bar is divided in 35 step indicators. For each required electrical test there is one step indicator provided.

- When an electrical test is finished tap the [Step done] button. Then the step indicator in the progress bar turns from a green frame to a green box.

The progress bar does not have any impact on the functionality of the electrosurgical generator. The values that are measured during the tests are not stored within the device.

All values and results of the safety check need to be recorded in the inspection card.

Relays are set automatically

For certain measurements closed output relays are needed. Closed output relays disable the electrosurgical generator for normal operation. During the assisted safety test, the electrosurgical generator automatically sets the relays as required for a certain test.

It is necessary to follow the preset sequence of steps and that the step number that is displayed on the touchscreen always corresponds with the test to be performed actually. Otherwise, the electrosurgical generator can be damaged seriously.

Other settings

The load resistance [Ohm] for the corresponding test must be set manually. The settings regarding socket, mode and power/effect are set automatically by the safety test assistant.

6.1.1 Starting the safety test assistant

Preconditions

- The electrosurgical generator is connected to the mains electricity.
- The electrosurgical generator is switched on.

1. Tap [Settings] on the home screen.
The [Settings] screen is displayed.
2. Tap [Maintenance].
The [Maintenance] tab is displayed.


3. Tap the [Safety test] button.
The first screen [Introduction] of the safety test assistant is displayed.
4. Consider the given information on the first screen.
5. If the progress bar displays green boxes, e.g. from the last year's test, then tap the [Reset all steps] button.
6. Tap the [Step done] button.
7. Tap the [Next] button to continue.

6.2 Electrical safety

The following 3 measurements belong to step 2 of the safety test assistant.

6.2.1 Checking the earth resistance and the earth leakage current

1. Switch off the electrosurgical generator.
2. Disconnect the power cable from the mains electricity.
3. Set up the test as shown in the figure above.
4. Ensure that power switch of the electrosurgical generator is depressed.
5. Verify that the protective earth resistance is $\leq 0.3 \Omega$.
6. Record the test results.
7. Verify under normal condition (NC) that the earth leakage current is $\leq 0.30 \text{ mA}$.
8. Verify under single fault condition (SFC) that the earth leakage current is $\leq 1.0 \text{ mA}$.
9. Record the test results.
10. At the electrical safety tester switch back to "normal condition" so that the ESG-410 turns on again.

6.2.2 Checking the patient leakage current for HF output

NOTICE


Serious damage due to wrong test step

Interconnecting all receptacles like it is needed for this measurement can seriously damage the electrosurgical generator.

Before interconnecting the 3 monopolar receptacles, ensure that the safety test assistant displays the test step 2/35.

Only for this measurement (step 2/35) it is allowed to interconnect the 3 receptacles of the MONOPOLAR socket.


1. Verify that the safety test assistant displays the appropriate step for this test.
2. Set up the test as shown in the figure above. Connect the cables to all 3 sockets of J059 or to J008.
3. Short circuit all connectors.
4. Ensure that the power switch of the electrosurgical generator is depressed. Consider that the device turns on only when the electrical safety tester is started.
5. Under normal condition (NC) for AC:
Verify that the patient leakage current is $\leq 10 \mu\text{A}$.
6. Under normal condition (NC) for DC:
Verify that the patient leakage current is $\leq 10 \mu\text{A}$.
7. Under single fault condition (SFC, "open earth") for AC:
Verify that the patient leakage current is $\leq 50 \mu\text{A}$.
8. Under single fault condition (SFC, "open earth") for DC:
Verify that the patient leakage current is $\leq 50 \mu\text{A}$.
9. Record the test results.
10. Disconnect all banana cables from the electrosurgical generator.
11. Tap [Step done] to confirm the measurements.
12. Tap the [Next] button to continue.

6.3 Checking the DC resistance – Step 3

Preconditions

- The electrosurgical generator is connected to the mains electricity.
- The electrosurgical generator is switched on.
- A compatible foot switch is connected.

1. Set up the test as shown in the figure above.
2. Verify that the resistance is $\geq 2 \text{ M}\Omega$.
3. Record the test results.
4. Dismantle the complete test setup.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.

6.4 Measurements for the HF leakage current

The following sections describe the measurements for the HF leakage current for each output socket of the electrocautery generator.

The HF current is measured as follows:


- With the electrocautery analyzer through an internal resistance of 200Ω .
- From the appropriate output connector of each pole.
- The electrocautery generator is operated at maximum output power settings in an appropriate mode.
- The output is at rated load.

The load resistor $R_L = 500 \Omega$ (J005) is connected via the connection cable (J007, banana) with the active pole MONOPOLAR and via the connection cable J011 with the NEUTRAL connector. Alternatively (if supported), the “external load” function of the electrocautery analyzer can be used.

6.5 Checking the HF leakage current for MONOPOLAR 1


6.5.1 Loaded condition

At the active electrode terminal – Step 4

1. Set up the test as shown in the figure above with a load resistance of $R_L = 500 \Omega$.
2. Press any of the foot switch pedals.
3. Verify that the HF leakage current is $\leq 150 \text{ mA}$.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.


At the neutral electrode terminal – Step 5

1. Set up the test as shown in the figure above with a load resistance of $R_L = 500 \Omega$.
2. Press any of the foot switch pedals.
3. Verify that the HF leakage current is $\leq 150 \text{ mA}$.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.

At the active electrode terminal – Step 6


This is an optional test.

1. Set up the test as shown in the figure above with a load resistance of $R_L = 500 \Omega$.
2. Press any of the foot switch pedals.
3. Verify that the HF leakage current is $\leq 150 \text{ mA}$.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.

At the neutral electrode terminal – Step 7

This is an optional test.


1. Set up the test as shown in the figure above with a load resistance of $R_L = 500 \Omega$.
2. Press any of the foot switch pedals.
3. Verify that the HF leakage current is $\leq 150 \text{ mA}$.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.

6.5.2 Unloaded condition

At the active electrode terminal – Step 8


This is an optional test.

1. Set up the test as shown in the figure above.
2. Press any of the foot switch pedals.
3. Verify that the HF leakage current is ≤ 150 mA.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.

At the neutral electrode terminal – Step 9

This is an optional test.


1. Set up the test as shown in the figure above.
2. Press any of the foot switch pedals.
3. Verify that the HF leakage current is ≤ 150 mA.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.

6.6 Checking the HF leakage current for MONOPOLAR 2


6.6.1 Loaded condition

At the active electrode terminal – Step 10

1. Set up the test as shown in the figure above with a load resistance of $R_L = 500 \Omega$.
2. Press any of the foot switch pedals.
3. Verify that the HF leakage current is $\leq 150 \text{ mA}$.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.


At the neutral electrode terminal – Step 11

1. Set up the test as shown in the figure above with a load resistance of $R_L = 500 \Omega$.
2. Press any of the foot switch pedals.
3. Verify that the HF leakage current is $\leq 150 \text{ mA}$.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.


6.6.2 Unloaded condition

At the active electrode terminal – Step 12

1. Set up the test as shown in the figure above.
2. Press any of the foot switch pedals.
3. Verify that the HF leakage current is ≤ 150 mA.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.

At the neutral electrode terminal – Step 13


1. Set up the test as shown in the figure above.
2. Press any of the foot switch pedals.
3. Verify that the HF leakage current is ≤ 150 mA.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.

6.7 Checking the HF leakage current for BIPOLAR


6.7.1 Loaded condition

At the active electrode terminal 1 (left) – Step 14

1. Set up the test as shown in the figure above with a load resistance of $R_L = 500 \Omega$.
2. Press any of the foot switch pedals.
3. Verify that the HF leakage current is $\leq 70.7 \text{ mA}$.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.

At the active electrode terminal 2 (right) – Step 15


1. Set up the test as shown in the figure above with a load resistance of $R_L = 500 \Omega$.
2. Press any of the foot switch pedals.
3. Verify that the HF leakage current is $\leq 70.7 \text{ mA}$.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.

6.8 Checking the HF leakage current for UNIVERSAL 1


6.8.1 Loaded condition

At the active electrode terminal 1 (black) – Step 16

1. Set up the test as shown in the figure above with a load resistance of $R_L = 75 \Omega$. Connect the cable to AE of J059 or to J008.
2. Press any of the foot switch pedals.
3. Verify that the HF leakage current is $\leq 70.7 \text{ mA}$.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.


At the active electrode terminal 2 (red) – Step 17

1. Set up the test as shown in the figure above with a load resistance of $R_L = 75 \Omega$. Connect the cable to NE of J059 or to J008.
2. Press any of the foot switch pedals.
3. Verify that the HF leakage current is $\leq 126.5 \text{ mA}$.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.


6.8.2 Unloaded condition

At the active electrode terminal 1 (black) – Step 18

1. Set up the test as shown in the figure above. Connect the cable to AE of J059 or to J008.
2. Press any of the foot switch pedals.
3. Verify that the HF leakage current is ≤ 126.5 mA.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.

At the active electrode terminal 2 (red) – Step 19


1. Set up the test as shown in the figure above. Connect the cable to NE of J059 or to J008.
2. Press any of the foot switch pedals.
3. Verify that the HF leakage current is ≤ 126.5 mA.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.

6.9 Checking the HF leakage current for UNIVERSAL 2


6.9.1 Unloaded condition - HF leakage current for UNIVERSAL 2

At the active electrode terminal 1 (black) – Step 20

1. Set up the test as shown in the figure above. Connect the cable to AE of J059 or to J008.
2. Press any of the foot switch pedals.
3. Verify that the HF leakage current is ≤ 100 mA.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.

At the active electrode terminal 2 (red) – Step 21

1. Set up the test as shown in the figure above. Connect the cable to NE of J059 or to J008.
2. Press any of the foot switch pedals.
3. Verify that the HF leakage current is ≤ 100 mA.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.
7. Disconnect the electrosurgical analyzer (J002) from the potential equalization connector.

6.10 Checking the hand switch functions – Step 22

For this test, the hand switch functions must be checked as follows:

- For the socket MONOPOLAR 1 with tool J009: cutting and coagulation
- For the socket MONOPOLAR 2 with tool J009: cutting and coagulation
- For the socket UNIVERSAL 1 with tool J059: cutting and coagulation
- For the socket UNIVERSAL 2 with tool J059: cutting and coagulation

1. Connect the corresponding tool to the output socket.
2. Press the hand switch.
3. Verify that the box [Press hand switch] on the touchscreen takes on the correct background color, i.e.:
 - yellow for cutting
 - blue for coagulation
4. Record the test result.
5. Continue with the next hand switch function.
6. When all hand switch functions are checked, tap [Step done] to confirm the test.
7. Tap the [Next] button to continue.

6.11 Checking the HF output power, open circuit and super ignition

NOTICE

If a 4 mm plug is connected incorrectly to the MONOPOLAR socket, this socket can be destroyed during activation.

Only connect a 4 mm plug to the receptacle on the right side.

For the measurements of these high current values, it is important, that the banana cables do not act as inductors.


Make sure that the cables are uncoiled completely.

6.11.1 Checking the HF output power for BIPOLAR – Steps 23-25

Power limits for BIPOLAR

Step	Mode	Power	Effect	Load	Requirement
23	AutoCoag	120	3	75 Ω	$96.8 \text{ W} \leq P_{\text{out}} \leq 145.2 \text{ W}$
24	HardCoag	120	3	500 Ω	$7.44 \text{ W} \leq P_{\text{out}} \leq 11.1 \text{ W}$
25	HardCoag	120	3	5000 Ω	Message i0001 "Open circuit" is displayed

Test description


1. Set up the test as shown in the figure above.
2. Press any of the foot switch pedals.
3. Verify that the output power is within the limit according to the table above.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.
7. Repeat this test for the safety test assistant steps 24 and 25.

6.11.2 Checking the HF output power for MONOPOLAR 1 – Steps 26-30

Power limits for MONOPOLAR 1

Step	Mode	Power	Effect	Load	Requirement
26	SoftCoag	200	5	50 Ω	$158.4\text{ W} \leq P_{\text{out}} \leq 237.6\text{ W}$
27	SprayCoag	120	3	500 Ω	$98.4\text{ W} \leq P_{\text{out}} \leq 147.6\text{ W}$
28	PowerCoag	120	3	500 Ω	$96.0\text{ W} \leq P_{\text{out}} \leq 144.0\text{ W}$
29	ForcedCoag	120	3	500 Ω	$94.4\text{ W} \leq P_{\text{out}} \leq 141.6\text{ W}$
30	ForcedCoag	120	3	2000 Ω	$28.0\text{ W} \leq P_{\text{out}} \leq 42.0\text{ W}$

Test description


1. Set up the test as shown in the figure above.
2. Press any of the foot switch pedals.
3. Verify that the output power is within the limit according to the table above.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.
7. Repeat this test for the safety test assistant steps 27 and 30.

6.11.3 Checking the HF output power for MONOPOLAR 2 – Steps 31-33

Power limits for MONOPOLAR 2

Step	Mode	Power	Effect	Load	Requirement
31	Blend	200	–	300 Ω	$159.2 \text{ W} \leq P_{\text{out}} \leq 238.8 \text{ W}$
32	Pure	300	–	300 Ω	$236.8 \text{ W} \leq P_{\text{out}} \leq 355.2 \text{ W}$
33	Fulgurate	120	–	500 Ω	$95.2 \text{ W} \leq P_{\text{out}} \leq 142.8 \text{ W}$

Test description


1. Set up the test as shown in the figure above.
2. Press any of the foot switch pedals.
3. Verify that the output power is within the limit according to the table above.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.
7. Repeat this test for the safety test assistant steps 32 and 33.

6.11.4 Checking the super ignition for UNIVERSAL 1 – Step 34

Power limits for UNIVERSAL 1

Step	Mode	Power	Effect	Load	Requirement
34	SuperIgnition test	100	1	0 Ω	$2.4\text{ A} \leq I_{out} \leq 3.5\text{ A}$ and no error occurs


1. Set up the test as shown in the figure above. Connect the cable to AE and NE of J059 or to J008.
2. Press any of the foot switch pedals.
3. Verify that the current ranges between 2.4 and 3.5 A.
4. Verify that no error occurs.
5. Record the test results.
6. Tap [Step done] to confirm the measurements.
7. Tap the [Next] button to continue.

6.11.5 Checking the HF output power for UNIVERSAL 1 – Step 35

Power limits for UNIVERSAL 1

Step	Mode	Power	Effect	Load	Requirement
35	PK PureCut	200	3	200 Ω	$153.6 \text{ W} \leq P_{\text{out}} \leq 230.4 \text{ W}$

Test description

1. Set up the test as shown in the figure above. Connect the cable to AE and NE of J059 or to J008.
2. Press any of the foot switch pedals.
3. Verify that the output power is within the limit according to the table above.
4. Record the test results.
5. Tap [Step done] to confirm the measurements.
6. Tap the [Next] button to continue.

End of assistance for electrical safety test

With this measurement the assistance for the electrical safety test ends. Also observe the following tests to complete the regular safety check.

7 Contact quality monitoring (CQM)

Preconditions

- The electrosurgical generator is connected to the mains electricity.
- The foot switch is connected to the electrosurgical generator.

7.1 CQM test for bipolar modes

1. Switch on the electrosurgical generator.
2. Ensure that the neutral electrode socket is left open, i.e. nothing is connected.
3. Assign the foot switch to the BIPOLAR socket.
4. Verify that the CQM indicator with the status "No neutral electrode connected" is displayed on the touchscreen.
5. Select any settings for any bipolar mode.
6. Press one of the foot switch pedals.
7. Verify that activation is possible.
8. Record the test result.

7.2 CQM test for monopolar modes

Simulating a non-split type neutral electrode

1. Connect the neutral electrode cable J011 to the [NEUTRAL] socket.
2. Short-circuit the neutral electrode cable J011, 0 Ω .
3. Verify that a message box regarding non-split neutral electrodes is displayed.
4. On the home screen, verify that the indicator for the non-split type neutral electrode is displayed in grey.
5. Connect the load resistor J004 with a load of 5 Ω to the neutral electrode socket.
6. Verify that a message box regarding non-split neutral electrode is displayed.
7. On the home screen, verify that the indicator for the non-split type neutral electrode is displayed in grey.
8. Record the test result.

Simulating a split type neutral electrode

9. Disconnect the neutral electrode cable J011.
10. Change the load of the load resistor J004 to 140 Ω .
11. Reconnect the neutral electrode cable J011 to the [NEUTRAL] socket.
12. On the home screen, verify that the indicator for the split type neutral electrode is displayed in green.
13. Record the test result.

8 Final tests

8.1 Checking the power switch, sound volume and display

LED of power switch

1. Verify that the circular LED of the power switch shines green when the electrosurgical generator is switched on.

Volume

2. Tap [Settings] on the home screen.
The [Setup] tab of the [Settings] screen is displayed.
3. In the volume pane tap [+] to increase the volume. Tap [-] to decrease the volume.
4. Verify that the volume of the tone increases and decreases accordingly.
5. Adjust the volume by the volume control at the rear panel of the electrosurgical generator.
6. Verify that the volume of the tone increases and decreases accordingly.

Brightness

7. In the brightness pane tap [+] to increase the brightness. Tap [-] to decrease the brightness.
8. Verify that the brightness increases and decreases accordingly.
9. Record the test results.

8.2 Enhanced system check

The enhanced system check tests if the electrosurgical generator is able to correctly recognize the connections at its sockets.

8.2.1 Testing without connections

As a first step, the test is done in the condition that nothing is connected to any of the sockets:

1. Check that no test tool is connected to any of the socket.
2. Tap [Settings] on the home screen.
The [Settings] screen is displayed.
3. Tap the [Maintenance] tab.
The [Maintenance] tab is displayed.
4. Tap the [System check] button.
The message box [Enhanced system check] is displayed.
5. Tap [OK] to start the system test.
The electrosurgical generator starts the enhanced system check.
6. Confirm that the message box [System check completed successfully] is displayed and that for each socket the statement is [Not connected.]
7. Record the test result.

8.2.2 Testing with connections

After the test without connections, the enhanced system check is done for each of the sockets of the electrosurgical generator.

Connections for the foot switch sockets

For this test 2 different foot switch types are required, i.e.:

- 1 single-pedal foot switch (wired) and 1 double- pedal foot switch (wired) or
- 1 single-pedal foot switch (wired) and 1 double-pedal foot switch (wireless) or
- 1 double-pedal foot switch (wired) and 1 double-pedal foot switch (wireless)

The test must be performed with each foot switch type at each of the 3 foot switch sockets.

Connections for the sockets on the front panel

- UNIVERSAL 1 socket with J008/J059
- UNIVERSAL 2 socket with J008/J059
- Neutral electrode socket with J002 (short circuit)

Performing the test

The test with connections can be done for all sockets at the same time or they can be done one after another for each socket separately.

1. Check that a test tool is connected to the corresponding socket.
2. Tap [Settings] on the home screen.
The [Settings] screen is displayed.
3. Tap the [Maintenance] tab.
The [Maintenance] tab is displayed.
4. Tap the [System check] button.
The message box [Enhanced system check] is displayed.
5. Tap [OK] to start the system test.
The electrosurgical generator starts the enhanced system check.
6. Confirm that the message box [System check completed successfully] is displayed and that the statement for each socket matches the connected test tool.
7. Record the test result.

8.3 Foot switch function test

For this test a double-pedal foot switch is required. The test must be performed for all 3 foot switch sockets.

1. Connect a double-pedal foot switch to the first foot switch socket.
2. Assign the connected double-pedal foot switch to the BIPOLAR socket.
3. Confirm that the corresponding foot switch indicator is displayed in the bipolar pane on the home screen.
4. Press the blue coag pedal.
5. Confirm that the background color of the bipolar pane changes to blue.
6. Press the yellow cut pedal.
7. Confirm that the background color of the bipolar pane changes to yellow.
8. Repeat this test for the other 2 foot switch sockets.
9. Record the test result.

8.4 Communication test

1. Establish a ring communication by connecting the communication cable (J012) to the LINK-IN socket and the LINK-OUT socket.
2. Verify that the error message [E0394 Communication error] is displayed on the touchscreen.
3. Record the test result.

8.5 USB socket test

1. Connect an USB flash drive (J021) to the USB socket.
2. Tap [Settings] on the home screen.
The [Settings] screen is displayed.
3. Tap the [Maintenance] tab.
The [Maintenance] tab is displayed.
4. Tap the [Data transfer] button.
The [Data transfer] screen is displayed.
5. Tap the [Export error log] button.
The progress of the data export is displayed.
6. Verify that the export is completed successfully.
7. Record the test result.

8.6 Date and time


1. Tap [Settings] on the home screen.
The [Settings] screen is displayed.
2. For [Date & time & format] check if the date and time is displayed correctly.

If the date and time need correction, proceed as follow:

1. Tap the options button for [Date & time & format].
The [Date & time & format] screen is displayed.
2. Adjust the settings.
3. Tap [OK] to save the settings.
The display returns to [Settings] screen. The updated date/time/format settings are displayed.

8.7 Pairing status of wireless foot switch


This test is applicable only for wireless foot switches.

1. Check that the pairing stickers on the foot switch and on the receiver match, as shown in the sample figure above.
2. With the serial cable, connect the receiver to one of the foot switch sockets of the electro-surgical generator.
3. Switch on the electro-surgical generator.
4. Press a foot switch pedal.
5. Ensure that the connection indicator on the receiver shines green.
6. Record the test result.

8.8 Inspection label

Only applicable for certain countries.

1. Stick a new inspection label each to the rear panel of the electrosurgical generator and to the rear of the foot switch.
2. Mark the inspection labels with the due dates for the next safety check.
3. Record the due dates in the inspection card.

9 Disconnecting

- Disconnect the devices from the mains electricity.
- Disconnect all test tools.

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WA91301W
WA91303W

Contents

1 General information	5
1.1 User instructions	5
1.2 Conventions throughout this document.....	5
1.3 Safety-related information in this document.....	6
2 Safety information.....	8
2.1 Intended purpose	8
2.2 Contact.....	9
2.3 General safety notes	10
3 Product description	11
3.1 Scope of delivery	11
3.2 Operating elements.....	12
3.3 Symbols	14
3.4 Warranty.....	16
4 Installation	17
4.1 Connecting the foot switch	17
4.2 Placing the foot switch	18
5 Before use.....	19
5.1 Inspection.....	19
6 Use.....	21
6.1 Adapting the foot switch position.....	21
6.2 Activation of the electrosurgical generator	21
6.3 Changing the foot switch assignment	21
7 After use	22
7.1 Disconnecting	22
8 Reprocessing	23
8.1 General information for reprocessing.....	23
8.2 Cleaning	23
8.3 Disinfection.....	24

9 Storage and maintenance	25
9.1 Storage.....	25
9.2 Periodic safety checks	25
10 Repair, shipment and disposal	26
10.1 Repair	26
10.2 Shipment.....	27
10.3 Disposal	27
11 Compatible equipment	28
12 Technical data	29
12.1 Specifications for WA91301W (double pedal).....	29
12.2 Specifications for WA91303W (single pedal)	29
12.3 Ambient conditions.....	29
13 Index.....	30

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Element	Explanation
	This is the safety alert symbol. It is used to alert the user to potential physical injury hazards. Observe all safety messages that follow this symbol to avoid possible injury.
	This symbol indicates additional helpful information.
1. 2. 3.	A numeration indicates a sequence of actions.
●	Bullet points indicate individual actions.
▶	Arrows indicate the results of actions.
— — —	Dashes indicate the listing of data, options or objects.

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Repair and product modifications

Unauthorized repairs and product modification can cause injury to the patient and the user. Repairs may only be performed by qualified servicing personnel that has been authorized by Olympus. Refer to the section "10 Repair, shipment and disposal" on page 26 for more details.

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
Contents


- 1) Foot switch (single pedal or double pedal)
- 2) Instructions for use and contact information for customers

3.2 Operating elements


3.2.1 Single-pedal foot switch

Top view

- 1) [COAG] pedal, blue
To activate the selected coagulation mode.
- 2) Connection cable with foot switch plug


Bottom view

- 1) Identification plate


3.2.2 Double-pedal foot switch

Top view

- 1) Connection cable with foot switch plug
- 2) Toggle button
To change the foot switch assignment remotely.
- 3) [COAG] pedal, blue
To activate the selected coagulation mode.
- 4) Bail handle
Engages at 0°, 90°, 180° and 270°.
- 5) [CUT] pedal, yellow
To activate the selected cutting mode.

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
Symbol	Explanation
MD	Medical device
SN	Serial number
QTY	Quantity
REF	Catalog number
Rx only	Federal (USA) law restricts this device to sale by or on the order of a physician
CE	CE marking indicating that the device is in conformity with the applicable requirements set out in applicable European Union harmonization legislation

Symbol	Explanation
	Manufacturer
	Caution, consult accompanying documents
	Consult instructions for use
	Unique device identifier
	Transport conditions
	Storage conditions
	Temperature limit
	Humidity limit
	In accordance with European Directive on Waste Electrical and Electronic Equipment (WEEE), this symbol indicates that the product must not be disposed of as unsorted municipal waste, but should be collected separately
	General symbol for recovery/recyclable
	Green Dot symbol for dual recycling system
IPX8	Protection class according to IEC 60529: Protected against the effects of continuous immersion in water

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4 Installation

This chapter explains how to install the foot switch. This includes the information on how to place the foot switch correctly.

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
4.1 Connecting the foot switch

1. Correctly align the foot switch plug with the foot switch socket on the rear panel of the electrosurgical generator. Observe the red markings.
2. Insert the plug into the foot switch socket until the plug clicks into place.

- ▶ The connection is fixed.

4.2 Placing the foot switch

- Place the foot switch on an even surface.
- Do not place the foot switch in corridors or in areas that people pass during the procedure.

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5.1.1 Inspecting for visible damage

1. Check that the foot switch has:
 - no dents, cracks, kinks, or deformations
 - no deep scratches
 - no corrosion
 - no missing or loose parts
2. Check that the foot switch and the receiver are visually clean.
3. Check all markings on the foot switch and the receiver for clear visibility.
4. Check that the pedals and the toggle button of the foot switch run smoothly.
5. Do not use the foot switch if any kind of damage or problem is observed.

5.1.2 Inspecting for functionality

1. Switch on the electrosurgical generator.
2. Assign the foot switch to an output socket of the electrosurgical generator. Refer to the instructions for use of the electrosurgical generator.

3. Ensure that the corresponding foot switch indicator is displayed on the touchscreen.
4. Press one pedal of the foot switch.
5. Check that the activation tone of the electrosurgical generator sounds.
6. Do not use the foot switch if any kind of problem is observed.

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6.2 Activation of the electrosurgical generator

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6.3 Changing the foot switch assignment

Use the toggle button to change the foot switch's assignment to an output socket of the electrosurgical generator. With each tap on the toggle button the assignment switches clockwise to the next output socket.

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7.1 Disconnecting

1. Switch off the electrosurgical generator.
2. At the foot switch plug, pull back the outer sleeve of the connector and demate the plug.

8 Reprocessing

8.1 General information for reprocessing

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1. Take a soft, lint-free cloth that is soaked with a suitable reprocessing agent.

2. Thoroughly wipe the surfaces of the foot switch until visibly clean.
3. In case of heavy soiling, use additional cloths.
4. Before disinfection, let the surfaces dry or wipe them dry using a lint-free cloth.

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Disinfection procedure

1. Take a soft, lint-free cloth that is soaked with a suitable disinfectant.
2. Thoroughly wet the surfaces of the foot switch by wiping with moderate pressure.
3. In a wet state, completely depress each of the pedals and wipe each pedal with moderate pressure
4. The surfaces must remain visibly wet for the contact time that is designated by the disinfectant's manufacturer but not less than 3 minutes.
5. Use additional cloths to assure the designated continuous wet contact time.
6. Let the surfaces dry or wipe them dry using a lint-free cloth.

9 Storage and maintenance

9.1 Storage

The foot switch must be stored according to the environmental conditions as described in the section "12.3 Ambient conditions" on page 29.

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
System chart

12 Technical data

12.1 Specifications for WA91301W (double pedal)

Dimensions	392 x 231 x 49 mm
Weight	1,600 g
Voltage	max. AC 25 V/DC 60 V
Cable length	4 m
Degree of protection (IEC 60529)	IPX8 (n/a for plug, cable)
Expected service life	10 years

12.2 Specifications for WA91303W (single pedal)

Dimensions	152 x 161 x 38 mm
Weight	830 g
Voltage	max. AC 25 V/DC 60 V
Cable length	4 m
Degree of protection (IEC 60529)	IPX8 (n/a for plug, cable)
Expected service life	10 years

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Storage conditions

Temperature	10 to 40 °C
Relative humidity	30 to 85%
Atmospheric pressure	50 to 160 hPa (≤ 3000 m)

Transport conditions

Temperature	-40 to 70 °C
Relative humidity	10 to 90%
Atmospheric pressure	50 to 160 hPa (≤ 3000 m)

13 Index

A	
Annual safety checks	25
C	
CAUTION, signal word	7
Cleaning agents	23
Clinical benefits	9
Contraindications.....	8
D	
DANGER, signal word	7
Disinfectants.....	24
G	
General safety notes	6
I	
Incident reporting	9
Intended use.....	8
Intended user	8
i-symbol.....	5
M	
Manufacturer	9
N	
NOTICE, signal word.....	7
O	
Operating conditions	29
P	
Patient target group.....	8
R	
Reprocessing agents.....	23

S

Safety alert symbol.....	5
Safety checks	25
Safety messages	7
Service centers	26
Signal word.....	7
Specific safety notes	6
Storage conditions.....	29
System chart	28

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CE₀₁₉₇ WA91302W

Contents

1 General information	5
1.1 User instructions	5
1.2 Conventions throughout this document.....	5
1.3 Safety-related information in this document.....	6
2 Safety information.....	8
2.1 Intended purpose	8
2.2 Contact.....	9
2.3 General safety notes	10
3 Product description	11
3.1 Scope of delivery	11
3.2 Operating elements.....	12
3.3 LED indicators.....	14
3.4 Symbols	17
3.5 Warranty.....	19
4 Unpacking and distinctive marking.....	20
4.1 Unpacking	20
4.2 Attaching the pairing stickers	20
5 Installation	23
5.1 Safety notes for installation	23
5.2 Insertion of batteries	23
5.3 Connection of the receiver	26
5.4 Placement of foot switch and receiver	27
6 Before use.....	29
6.1 Safety notes for before use	29
6.2 Inspection.....	29
6.3 Consider other radio systems	31
7 Use.....	32
7.1 Adapting the foot switch position.....	32
7.2 Activation of the electrosurgical generator	32

7.3 Changing the foot switch assignment	33
8 After use	34
8.1 Disconnecting the receiver	34
9 Reprocessing	35
9.1 General information for reprocessing	35
9.2 Cleaning	35
9.3 Disinfection.....	36
10 Storage and maintenance	37
10.1 Storage.....	37
10.2 Periodic safety checks	37
10.3 Changing the batteries	38
11 Repair, shipment and disposal	39
11.1 Repair	39
11.2 Shipment.....	40
11.3 Disposal	40
12 Compatible equipment	41
13 Technical data	42
13.1 Specifications for the foot switch.....	42
13.2 Specifications for the receiver	42
13.3 Radio specifications	42
13.4 Ambient conditions.....	42
14 Index	44

1 General information

This section provides general information on how to use these instructions for use.

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Element	Explanation
	This is the safety alert symbol. It is used to alert the user to potential physical injury hazards. Observe all safety messages that follow this symbol to avoid possible injury.
	This symbol indicates additional helpful information.
1. 2. 3.	A numeration indicates a sequence of actions.
●	Bullet points indicate individual actions.
►	Arrows indicate the results of actions.
— — —	Dashes indicate the listing of data, options or objects.

Element	Explanation
1) 2) 3)	Numbers with right parenthesis name elements in illustrations

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2.1 Intended purpose

2.1.1 Intended use

● Foot switch to activate and deactivate the power output of a compatible electrosurgical generator.
Do not use for any other purpose.

2.1.2 Contraindications

There are no known contraindications. The intended use and contraindications of the equipment used in combination with the foot switch must be observed.

2.1.3 Patient target group

● Not dedicated to any patient population.

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Repair

Repair of the product may only be performed by trained qualified servicing personnel that has been authorized by Olympus. Otherwise, Olympus cannot be held responsible for the safety and performance of the product.

2.1.5 Clinical benefits

The device is used as part of a system. The clinical benefit of the device lies in its functionality as described by its intended purpose.

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2.2.2 Incident reporting

European Union

If a serious incident occurs with the device, report it to the manufacturer and the relevant national authority.

Eurasian Economic Union

If an undesirable event occurs with the device, report it to the manufacturer or the authorized representative.

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2.3 General safety notes

The following safety information addresses points that apply to the general handling and use of the product. This information is to be supplemented by the information given in each chapter of this document.

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Repair and product modifications

Unauthorized repairs and product modification can cause injury to the patient and the user. Repairs may only be performed by qualified servicing personnel that has been authorized by Olympus. Refer to the section “11.1 Repair” on page 39 for more details.

Regular safety check

To ensure electrical safety and to maintain device performance, perform a safety check on all medical electrical equipment annually or according to national regulations and guidelines.

3 Product description


This section provides a general overview of the product and its functions. This includes the following information:

- Scope of delivery
- Operating elements
- Product functions
- Explanation of symbols
- Warranty information

3.1 Scope of delivery

Before use, check that all items listed below are available. Contact an Olympus representative or an authorized service center if any items are missing or damaged.

Contents

- 1) Wireless foot switch, double pedal
- 2) Receiver
- 3) 2 sheets with large pairing stickers (numbers 1-5)

- 4) 2 sheets with small pairing stickers (numbers 1-5)
- 5) Instructions for use and contact information for customers
- 6) Serial cable


Compatible batteries are not included in the scope of delivery. Refer to the section "5.2 Insertion of batteries" on page 23.

3.2 Operating elements


3.2.1 Double-pedal foot switch

Top view

- 1) Toggle button
To change the foot switch assignment remotely.
- 2) Battery indicator
- 3) Connection indicator
- 4) [COAG] pedal, blue
To activate the selected coagulation mode.
- 5) Bail handle
Engages at 0°, 90°, 180° and 270°.
- 6) [CUT] pedal, yellow
To activate the selected cutting mode.


Bottom view

- 1) Screw cap of battery compartment
- 2) Battery illustration
- 3) Identification plate

3.2.2 Receiver and serial cable


Receiver

- 1) Connection indicator
- 2) Power indicator

- 3) Serial port
To connect the serial cable
- 4) Battery indicator
- 5) Pairing button
For service purposes only. Do not depress.
- 6) Feet

Serial cable

- 1) Foot switch plug with push-pull locking
To connect the serial cable to the foot switch socket of the electrosurgical generator.
- 2) Receiver plug with bayonet locking
To connect the serial cable to the serial port of the receiver.

3.3 LED indicators

This section describes the meanings of the different LED indicator conditions.

3.3.1 Conditions at the foot switch

Connection indicator	Battery indicator	Description
Off	Off	Sleep mode, foot switch has not been activated in the last 3 seconds.

Connection indicator 	Battery indicator 	Description
Off	On	No pairing partner found. Possible reasons: The receiver is out of range. The receiver is not connected properly.
Green	(irrelevant)	Radio connection between foot switch and the receiver is established.
Green, blinking slowly	(irrelevant)	Radio connection between foot switch and the receiver is established, with limited connection quality.
Green, quickly blinking	(irrelevant)	For service purposes only.
Dark orange, blinking slowly	(irrelevant)	For service purposes only.
Green with dark orange pulses	(irrelevant)	For service purposes only.
Dark orange, blinking quickly	(irrelevant)	For service purposes only.
Dark orange	(irrelevant)	System error. Contact an authorized service center.
(irrelevant)	Green	Battery capacity $\geq 50\%$
(irrelevant)	Orange	Battery capacity between 20 and 50%.
(irrelevant)	Orange, blinking slowly	Battery capacity $< 20\%$. It is safe to finish the current procedure. Before starting a new procedure, change the batteries.
(irrelevant)	Orange, blinking quickly	Battery capacity below limit. The foot switch cannot be used.

3.3.2 Conditions at the receiver

Power indicator	Description
Green	Power supply is on.
Dark orange	System error. Contact an authorized service center



Connection indicator 	Battery indicator 	Description
Off	Off	Sleep mode, foot switch has not been activated for more than 3 seconds.
Green	(irrelevant)	Radio connection is established. Foot switch pedal is pressed.
Orange	(irrelevant)	Radio connection is established, no foot switch pedal is pressed.
Green with dark orange pulses	(irrelevant)	For service purposes only.
Dark orange, blinking slowly	(irrelevant)	For service purposes only.
Dark orange, blinking quickly	(irrelevant)	For service purposes only.
Dark orange	Off	Receiver is not paired with a foot switch. Contact an authorized service center.
Dark orange	Orange	System error. Contact an authorized service center.
(irrelevant)	Green	Battery capacity $\geq 50\%$
(irrelevant)	Orange	Battery capacity between 20 and 50%.
(irrelevant)	Orange, blinking slowly	Battery capacity $< 20\%$. It is safe to finish the current procedure. Before starting a new procedure, change the batteries.

Connection indicator 	Battery indicator 	Description
(irrelevant)	Orange, blinking quickly	Battery capacity below limit. The foot switch cannot be used.
Off, even is pedal is pressed	Orange	System error. Contact an authorized service center.
Dark orange	Orange	System error. Contact an authorized service center

3.4 Symbols

This section gives an explanation for each symbol used on the product and on the packaging of the product.

Symbol	Explanation
	Medical device
	Serial number
	Quantity
	Catalog number
Rx only	Federal (USA) law restricts this device to sale by or on the order of a physician
	CE marking indicating that the device is in conformity with the applicable requirements set out in applicable European Union harmonization legislation
	Manufacturer

Symbol	Explanation
	Refer to instructions for use
	Caution, consult accompanying documents
	Unique device identifier
	Transport conditions
	Storage conditions
	Temperature limit
	Humidity limit
	In accordance with European Directive on Waste Electrical and Electronic Equipment (WEEE), this symbol indicates that the product must not be disposed of as unsorted municipal waste, but should be collected separately
	General symbol for recovery/recyclable
	Green Dot symbol for dual recycling system
	Connection indicator
	Battery indicator
IPX8	Protection class according to IEC 60529: Protected against the effects of continuous immersion in water

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4 Unpacking and distinctive marking

Though not physically joined, a foot switch and its receiver are distinctively paired with each other by a unique software key that has been exchanged between the 2 devices. This is why a certain foot switch works only with a certain receiver.

When having more than one wireless foot switch, it easily happens that the several devices get mixed up. To avoid this, the devices must be distinctively marked with the provided pairing stickers.


The foot switch and its receiver are delivered in a paired condition. So there is no need for the user to take care of pairing.

4.1 Unpacking

When unpacking the foot switch and its receiver, ensure that the devices are not mixed up with other wireless foot switches.

- Only unpack one wireless foot switch at a time.
- Complete the marking with the pairing stickers before unpacking and marking another wireless foot switch.

4.2 Attaching the pairing stickers


- 1) Large pairing stickers for the foot switch
- 2) Small pairing stickers for the receiver

Olympus recommends to use two-digit numbers for the marking.

Large stickers for the foot switch

1. Peel off 2 of the large stickers, e.g. the numbers 1 and 1.


2. Attach the 2 large stickers to the top left of the foot switch so that they can be recognized easily.

Do not attach the stickers to the lower part of the foot switch, as they can abrade easily during use.

Small sticker for the receiver

1. Peel off the same 2 numbers from the sheet with the small stickers, e.g. 1 and 1.
2. Attach the 2 small stickers to the top of the receiver so that they can be recognized easily.

- Do not attach the stickers over the groove of the housing, as dust or cleaning liquids can cause the stickers to detach easily.
- Do not attach the stickers at the bottom of the receiver, as they will be hidden by the receiver holder of the cart.

Marking another foot switch and its receiver

If another wireless foot switch is setup, ensure that the foot switch and its receiver are also distinctively marked with different pairing stickers, e.g. 2 and 2.

5 Installation

This chapter explains how to install the foot switch. This includes the following information:

- Inserting the batteries
- Placing the foot switch correctly

5.1 Safety notes for installation

Precondition for installation

To avoid accidental activation, the electrosurgical generator must be switched off.

No rechargeable batteries

Only use single-use batteries for the foot switch. When mixing single-use batteries with rechargeable batteries, then the risk of a battery fire rises.

The discharge characteristics of rechargeable batteries are different from the discharge characteristics of single-use batteries. The different states of the LED indicators, as described in section “3.3 LED indicators” on page 14, are calculated based on single-use batteries. If rechargeable batteries are used, then it is possible that the LED indicator conditions regarding battery capacity do not work properly.

Mechanical stress

Do not expose the plug-connections to mechanical stress, as this can damage the electrical connections.

5.2 Insertion of batteries

Required batteries

3 pieces of

IEC no. code: R14

Nominal Ø: 26.2 mm

Nominal height: 50.0 mm

Common name: C or Baby


Olympus recommends to use one of the following battery types:

- Energizer Ultra+ LR14
- Panasonic Industrial LR14


Inserting the batteries

Olympus recommends to insert only completely new batteries.


1. Observe the battery illustration on the bottom of the foot switch.


2. Open the screw cap of the battery compartment on the rear of the foot switch.

3. Insert the 3 batteries considering the positive [+] and negative [-] poles.


4. Attach the screw cap to the battery compartment.

5. Close the screw cap of the battery compartment.


5.3 Connection of the receiver

Connecting the serial cable to the electrosurgical generator

1. Correctly align the foot switch plug with the foot switch socket on the rear panel of the electrosurgical generator. Observe the red markings.
2. Insert the plug into the foot switch socket until the plug clicks into place.
 - The connection is fixed.

Connecting the serial cable to the receiver

1. Correctly align the pins of the receiver plug with the serial port of the receiver. Observe the small arrow on the plug.
2. Insert the plug.
3. Turn the bayonet locking of the plug clockwise until it clicks into place.
 - The connection is fixed.
 - All 3 LEDs on the receiver light up orange for approximately 2 seconds. The foot switch is ready for use.
4. For correct placement of the receiver refer to the section “5.4 Placement of foot switch and receiver” on page 27.

The receiver is powered by the electrosurgical generator.

5.4 Placement of foot switch and receiver

Placing the receiver

- Place and operate the receiver in the non-sterile area only.
- Place the receiver so that the LEDs are visible during the procedure. E.g., the receiver can be placed on its feet on top of the electrosurgical generator.
- If available, mount the receiver on the receiver holder of the energy cart TC-E410. Refer to the instructions for use of the energy cart TC-E410.

If 2 wireless foot switches are connected to the electrosurgical generator:

- Ensure that the distance between the 2 receivers is at least 30 cm.
- If the receivers are mounted on the energy cart TC-E410, ensure that the 2 receiver holders are placed on opposite sides of the energy cart.

Placing the foot switch

- Place the foot switch on an even surface.
- Do not place the foot switch in corridors or in areas that people pass during the procedure.
- Place the foot switch within the maximum transmission range of the receiver, i.e. 10 m.

Transmission range

The following factors can reduce the transmission range:

- High-frequency interference
- WLAN networks or Bluetooth networks
- Mobile communication equipment, e.g. mobile phones
- Other wireless devices with the same frequency
- Inappropriate antenna alignment
- Metal objects
- Wall insulation coated with or containing metal
- Humans within the transmission range

Position sensor in the foot switch

As a safety feature, the foot switch cannot be activated when it is inclined by more than 35°. This prevents unintentional activation during pickup and storage of the foot switch.

6 Before use

This chapter describes all necessary steps to prepare the foot switch for each use.

6.1 Safety notes for before use

Sufficient battery capacity

The battery capacity is sufficient for the next procedure if the battery indicator on the foot switch shines green (capacity >50%) or orange without blinking (capacity between 20 and 50%).

If the battery indicator starts to blink orange slowly during a procedure, there is still sufficient capacity left to safely finish the current procedure.

If the battery indicator is orange and blinking quickly, then the foot switch cannot be used. The batteries must be changed.

6.2 Inspection

A thorough inspection ensures that the product is fit for use. Damaged products can cause an electric shock and thermal injury.

Replace if necessary

If the instrument is damaged or does not function properly, replace it. Contact an authorized service center for repair.

6.2.1 Inspecting for visible damage

1. Check that the foot switch has:
 - no dents, cracks, kinks, or deformations
 - no deep scratches
 - no corrosion
 - no missing or loose parts
2. Check that the foot switch and the receiver are visually clean.

3. Check all markings on the foot switch and the receiver for clear visibility.
4. Check that the pedals and the toggle button of the foot switch run smoothly.

Do not use the foot switch if any kind of damage or problem is observed.

6.2.2 Inspecting for functionality

Precondition

To avoid accidental activation, the electrosurgical generator must be switched off.

1. Check that the pairing stickers on the foot switch and on the receiver match.
 2. If the foot switch is in the sleep mode, press a pedal or the toggle button.
 3. Check that the battery capacity is higher than 20%, i.e. the battery indicators on the foot switch and on the receiver shine green or orange without blinking.
 - If the battery indicators are blinking orange, change the batteries before starting a new procedure.
 4. Ensure that the indicators of the receiver are visible during use.
 5. Switch on the electrosurgical generator.
 6. Assign the foot switch to an output socket of the electrosurgical generator. Refer to the instructions for use of the electrosurgical generator.
 7. Ensure that the corresponding foot switch indicator is displayed on the touchscreen.
 8. Press one pedal of the foot switch.
 9. Check that the activation tone of the electrosurgical generator sounds.
 10. Disconnect the receiver from the electrosurgical generator.
 11. Check that the foot switch indicator is not displayed anymore.
 12. Reconnect the receiver.
 13. Check that the foot switch indicator is displayed again.
- Do not use the foot switch if any kind of problem is observed.

6.3 Consider other radio systems

Other foot switches

When other wireless foot switches of the same frequency (2.4 GHz) are used within the same area and these foot switches transmit signals at the same time, transmission errors can occur. Do not use more than 6 wireless foot switches of this kind in the same area.

Bluetooth and WLAN

It is possible to operate the wireless foot switch in areas with Bluetooth and WLAN networks. If disturbances occur repeatedly, perform radio network planning.

7 Use

7.1 Adapting the foot switch position

If the bail handle is engaged in the upright position (90°), it can be used to conveniently move the foot switch by foot to the user's preferred position.

To minimize the risk of tripping, move the bail handle to its down position (0°) when the foot switch is not in use.

7.2 Activation of the electrosurgical generator

CAUTION

Risk of injury due to using the incorrect foot switch

When more than one electrosurgical generator is used with wireless foot switches, then activating the incorrect foot switch can lead to unintended activation of an incorrect electrosurgical generator.

- Before use, ensure to use the correct wireless foot switch. Check for matching pairing stickers on the foot switch and on the receiver.
- Make sure that the wireless foot switch is always near its receiver.
- Before the procedure, observe the instructions in the section “6.2.2 Inspecting for functionality” on page 30.

If the foot switch is not activated for more than 3 seconds, then it enters the sleep mode. The foot switch enters the active mode as soon as a pedal or the toggle button is depressed.

Activating cutting modes

All monopolar and bipolar cutting modes are activated by pressing the yellow [CUT] pedal of the foot switch.

Activating coagulation modes

All monopolar and bipolar coagulation modes are activated by pressing the blue [COAG] pedal of the foot switch.

Stopping activation

The HF output of the electrosurgical generator stops as soon as the foot switch pedal is released.

7.3 Changing the foot switch assignment

Use the toggle button to change the foot switch's assignment to another output socket of the electrosurgical generator. With each tap on the toggle button the assignment switches clockwise to the next output socket.

8 After use

This chapter describes all necessary steps that are required to prepare the foot switch for cleaning and disinfection.

8.1 Disconnecting the receiver

1. Switch off the electrosurgical generator.
2. At the foot switch plug, pull back the outer sleeve of the connector and demate the plug.
3. At the receiver plug, turn the bayonet locking anti-clockwise and demate the plug.

9 Reprocessing

Preconditions

- The electrosurgical generator must be switched off to avoid accidental activation.
- The battery compartment must be tightly closed to avoid ingress of liquids.

9.1 General information for reprocessing

The reprocessing cycle

- Reprocess the product according to the instructions in this chapter.
- Reprocess the product before each use.

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Cleaning procedure

Perform the following cleaning procedure after use. If necessary, let the product cool down. If cleaning is delayed, debris will solidify and it may be difficult to effectively clean the surfaces.

1. Take a soft, lint-free cloth that is soaked with a suitable reprocessing agent.
2. Thoroughly wipe the surfaces of the foot switch until visibly clean.
3. In case of heavy soiling, use additional cloths.
4. Before disinfection, let the surfaces dry or wipe them dry using a lint-free cloth.

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Inspecting before disinfection

1. Visually inspect the foot switch thoroughly after cleaning. The foot switch must be visually clean.
2. If there are any signs of debris, repeat the cleaning process.

Disinfection procedure

1. Take a soft, lint-free cloth that is soaked with a suitable disinfectant.
2. Thoroughly wet the surfaces of the foot switch by wiping with moderate pressure.
3. In a wet state, completely depress each of the pedals and wipe each pedal with moderate pressure.
4. The surfaces must remain visibly wet for the contact time that is designated by the disinfectant's manufacturer but not less than 3 minutes.
5. Use additional cloths to assure the designated continuous wet contact time.
6. Let the surfaces dry or wipe them dry using a lint-free cloth.

10 Storage and maintenance

10.1 Storage

The foot switch and its receiver must be stored according to the environmental conditions as described in the section “13.4 Ambient conditions” on page 42.

The energy cart TC-E410 provides a special compartment for storing foot switches.

10.2 Periodic safety checks

To ensure electrical safety, perform a periodic safety check on the foot switch and its receiver.

Frequency of safety checks

Perform the safety check annually or in accordance with applicable national regulations and guidelines.

Qualification for the safety check

The safety checks must only be performed by a qualified electrician with sufficient experience in maintaining medical electrical devices.

Appropriate environment and equipment

The safety check must only be performed in a technical laboratory with appropriate equipment.

Instructions in the maintenance manual

The detailed instructions for the periodic safety check are described in the maintenance manual for the electrosurgical generator ESG-410.

10.3 Changing the batteries

- Before changing the batteries, make sure that the electrosurgical generator is switched off.

Refer to the section “5.2 Insertion of batteries” on page 23.

11 Repair, shipment and disposal

This chapter provides instructions on the following:

- How to contact Olympus for any repairs
- How to ship the product to Olympus for repair or warranty claims
- Disposal of the product

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Sending used and contaminated products to Olympus represents a risk of infection. Olympus is entitled to refuse soiled or contaminated products for reasons of safety.

Reprocess the product before sending in.

If reprocessing will damage the product completely, clean the product as thoroughly as possible and mark it accordingly.

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
12 Compatible equipment

Olympus recommends the use of equipment as listed in this chapter. If combinations are used that are not listed in this chapter, the user takes the full responsibility.

Future equipment can also be compatible. For more information, contact an Olympus representative.

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System chart

13 Technical data

13.1 Specifications for the foot switch

Dimensions	392 x 211 x 49 mm
Weight	1,470 g
Voltage	4.5 V
Batteries	3 pieces of R14 (C/Baby)
Degree of protection (IEC 60529)	IPX8
Expected service life	10 years

13.2 Specifications for the receiver

Dimensions	102 x 76 x 38 mm
Weight	125 g
Length of serial cable	1.2 m
Voltage	max. AC 25 V/DC 60 V
Degree of protection (IEC 60529)	IPX8 (n/a for plug, cable)
Expected service life	10 years

13.3 Radio specifications

Antenna type	internal
Transmission range	10 m
Frequency	2.4 to 2.48 GHz
Channel spacing	2 MHz

13.4 Ambient conditions

Operating conditions

Ambient temperature	10 to 40 °C
Relative humidity	30 to 85%
Atmospheric pressure	70 to 160 hPa (≤ 3000 m)

Storage conditions

Temperature 10 to 40 °C
Relative humidity 30 to 85%
Atmospheric pressure 50 to 160 hPa (≤ 3000 m)

Transport conditions

Temperature -40 to 70 °C
Relative humidity 10 to 90%
Atmospheric pressure 50 to 160 hPa (≤ 3000 m)

14 Index

A	
Annual safety checks	37
B	
Batteries	23
Battery capacity	29
Bluetooth	31
C	
● CAUTION, signal word	7
Cleaning agents	35
Clinical benefits	9
Contraindications.....	8
D	
DANGER, signal word	7
Disinfectants	36
G	
General safety notes	6
I	
● Incident reporting	9
Intended use.....	8
Intended user	8
i-symbol.....	5
M	
Maintenance manual	37
Manufacturer	9
N	
NOTICE, signal word.....	7
Number stickers	20

O	
Operating conditions	42
P	
Pairing stickers	20
Patient target group	8
Position sensor	28
R	
Receiver	13
Reprocessing agents	35
S	
Safety alert symbol	5
Safety checks	37
Safety messages	7
Service centers	39
Signal word	7
Specific safety notes	6
Stickers	20
Storage conditions	43
System chart	41
T	
Toggle button, function	33
Transmission range	28
Transport conditions	43
W	
WARNING, signal word	7
WLAN	31

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Height	241 mm ±2 mm
Depth	501 mm ±2 mm
Depth including connected cables	520 mm ±10 mm
Weight	16,4 kg ±10%

Power supply:

Characteristic	Value
Voltage range	AC 100 to 120 V/AC 220 to 240 V
Frequency	50 Hz/60 Hz
Maximum input power for AC 100 to 120 V	1150 VA
Maximum input power for AC 220 to 240 V	1500 VA
Fuse	2 x 250 V/T 10 AH

Output:

Characteristic	Value
High frequency functions	Monopolar and bipolar
High frequency at rated load	344 to 516 kHz
Maximum HF power monopolar	300 W at 500 Ω
Maximum HF power bipolar	320 W at 75 Ω
All modes	25% duty cycle (e.g. activated for 10 s, deactivated for 30 s)

Plug requirements:

For the MONOPOLAR output sockets

Type of plug	Specification
3-pin plugs	diameter 4 mm
1-pin plugs	diameter 8 mm
1-pin plugs	diameter 4 mm
Coaxial plugs	diameter 9/5 mm

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Type of plug	Specification
2-pin plugs	pin diameter 4 mm, 28.6 mm pin spacing
2-pin plugs	pin diameter 4 mm with 3rd coding pin
2-pin plugs	pin diameter 4 mm, pin spacing 22 mm
Coaxial plugs	pin diameter 4/8 mm

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For the neutral electrode socket

Type of plug	Specification
2-pin plugs	diameter 2.5 mm, 10 mm pin spacing

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Type of electrode	Resistance
For split type neutral electrodes	from 25 to 150 Ohm
For non-split type neutral electrodes	<25 Ohm

Specifications for cables:

Power cables

Specification/type of cable	Value
WA95621A (for many European countries)	Type E/F
WA95622A (for Canada, USA and other countries)	Type B
WA95623A (for United Kingdom and other countries)	Type G
Length	4,5 m ±0,1 m
Plug on device side	angled

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Specifications for the USB flash drive:

Specification	Value
Storage capacity	max. 32 GB
Formatting	FAT32 or FAT

ExFAT, NTFS and others are not supported.
 Compatible with USB 1.1
 Multiple partitions are not supported.
 Do not connect hard disks drives.

Technical characteristics of Two-pedal switch WA91301W

Length	392 mm ± 0.5 mm
Width	231 mm ± 0.5 mm
Height	49 mm ± 0.3 mm
Weight	1600 g ± 10%
Length of cable	4 m ± 0.1 m
Degree of protection (IEC 60529)	IPX8 (n/a for plug, cable)

Technical characteristics of Single-pedal switch WA91303W

Length	152 mm ± 0.5 mm
Width	161 mm ± 0.5 mm
Height	38 mm ± 0.3 mm
Weight	830 g ± 10%
Length of cable	4 m ± 0.1 m
Degree of protection (IEC 60529)	IPX8 (n/a for plug, cable)

Technical characteristics of Two-pedal switch, wireless WA91302W

Length	392 mm ± 0.5 mm
Width	211 mm ± 0.5 mm
Height	49 mm ± 0.3 mm
Weight	1470 g ± 10%
Voltage	4,5 V
Batteries	3 pieces of R14 (C/Baby)

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Technical characteristics of Two-pedal switch, wireless WA91302W reciever

Length	102 mm ± 0.3 mm
Width	76 mm ± 0.3 mm
Height	38 mm ± 0.3 mm
Weight	125 g ± 10%
Length of serial cable	1,2 m ± 0.1 m
Voltage	max. AC 25 V/DC 60 V
Degree of protection (IEC 60529)	IPX8 (n/a for plug, cable)

Radio specifications

Antenna type	Internal
Transmission range	10 m
Frequency	2.4 to 2.48 GHz
Channel spacing	2 MHz

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HPCS – High Power Cut Support

This feature optimizes the start of the cutting procedure for certain cutting modes. By applying high power to the tissue, HPCS supports immediate spark ignition and

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reduces the risk of mechanical cutting.

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Description of the bipolar cutting modes:

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Control characteristics	HPCS
Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.3 to 1.8
Range of power levels	1 to 100 W

		PureCut
Коагуляция	AutoCoag	Fulgurate
	BiSoftCoag	ForcedCoag
	HardCoag	PowerCoag
	PK Coag	SoftCoag
	PK SoftCoag	Spray
	SalineCoag	SprayCoag

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PlasmaCut

Effect levels	1 to 3
Maximum output voltage at effect level 1	425 V _p
Maximum output voltage at effect level 2	485 V _p
Maximum output voltage at effect level 3	645 V _p
Maximum output current	1.5 A
Rated load	500 Ω

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Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.3 to 1.8
Range of power levels	1 to 320 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	597 V _p
Maximum output voltage at effect level 2	625 V _p
Maximum output voltage at effect level 3	700 V _p
Maximum output current	12 A
Rated load	75 Ω

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Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.3 to 1.8
Maximum application time	60 s

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Range of power levels	1 to 200 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	440 V _p
Maximum output voltage at effect level 2	520 V _p
Maximum output voltage at effect level 3	590 V _p
Maximum output current	3 A
Rated load	200 Ω
Limited activation time	60 s

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Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.3 to 1.8
Maximum application time	60 s
Range of power levels	1 to 300 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	500 V _p
Maximum output voltage at effect level 2	500 V _p
Maximum output voltage at effect level 3	560 V _p
Maximum output current	3 A
Rated load	75 Ω
Limited activation time	60 s

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PK PureCut

Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.3 to 1.8
Maximum application time	60 s
Range of power levels	1 to 200 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	410 V _p
Maximum output voltage at effect level 2	440 V _p
Maximum output voltage at effect level 3	510 V _p
Maximum output current	1.5 A
Rated load	200 Ω

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Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.3 to 1.8
Maximum application time	60 s
Range of power levels	1 to 200 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	430 V _p
Maximum output voltage at effect level 2	520 V _p
Maximum output voltage at effect level 3	610 V _p
Maximum output current	1.5 A
Rated load	200 Ω

Description of the bipolar coagulation modes:

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Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.3 to 1.8
Range of power levels	1 to 120 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	78 V _p
Maximum output voltage at effect level 2	134 V _p
Maximum output voltage at effect level 3	220 V _p
Maximum output current	1.4 A
Rated load	75 Ω

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Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.3 to 1.8
Range of power levels	1 to 120 W (max. 50 W with Autostart ON)
Effect levels	1 to 3
Maximum output voltage at effect level 1	95 V _p
Maximum output voltage at effect level 2	145 V _p
Maximum output voltage at effect level 3	220 V _p
Maximum output current	1.4 A
Rated load	75 Ω

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Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.3 to 1.8
Range of power levels	1 to 120 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	220 В _p
Maximum output voltage at effect level 2	220 В _p
Maximum output voltage at effect level 3	220 В _p
Maximum output current	3.3 А
Rated load	25 Ω

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Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.3 to 1.8
Range of power levels	1 to 120 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	150 В _p
Maximum output voltage at effect level 2	200 В _p
Maximum output voltage at effect level 3	220 В _p

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Maximum output current	2.3 A
Rated load	75 Ω

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Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.3 to 1.8
Range of power levels	1 to 120 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	150 V _p
Maximum output voltage at effect level 2	170 V _p
Maximum output voltage at effect level 3	230 V _p
Maximum output current	2.3 A
Rated load	75 Ω

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Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.3 to 1.8
Range of power levels	1 to 200 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	175 V _p
Maximum output voltage at effect level 2	230 V _p
Maximum output voltage at effect level 3	245 V _p

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level 3	
Maximum output current	5 A
Rated load	75 Ω

Description of the monopolar cutting modes:

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Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.8 to 2.9
Range of power levels	1 to 200 W
Effect levels	1
Maximum output voltage at effect level 1	2200 V _p
Maximum output current	1.0 A
Rated load	300 Ω

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Waveform	Continuous sinusoidal, unmodulated
Crest factor	3 to 4.5
Range of power levels	1 to 200 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	935 V _p
Maximum output voltage at effect level 2	1140 V _p
Maximum output voltage at effect level 3	1400 V _p
Maximum output current	1.5 A

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Control characteristics	HPCS
Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.3 to 1.8
Range of power levels	1 to 120 W
Effect levels	1
Maximum output voltage at effect level 1	560 В _p
Maximum output current	1.2 А
Rated load	500 Ω

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Waveform	Continuous sinusoidal, modulated
Crest factor	3.5 to 6
Range of power levels	1 to 120 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	1224 В _p
Maximum output voltage at effect level 2	1519 В _p
Maximum output voltage at effect level 3	2000 В _p
Maximum output current	1.2 А
Rated load	500 Ω

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Control characteristics	HPCS
Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.3 to 1.8
Range of power levels	20 to 150 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	697 V _p
Maximum output voltage at effect level 2	706 V _p
Maximum output voltage at effect level 3	714 V _p
Maximum output current	1.5 A
Rated load	100 Ω

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Control characteristics	HPCS
Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.3 to 1.8
Range of power levels	20 to 150 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	671 V _p
Maximum output voltage at effect level 2	704 V _p
Maximum output voltage at effect level 3	714 V _p

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Maximum output current	1.5 A
Rated load	100 Ω

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Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.3 to 1.8
Range of power levels	1 to 300 W
Effect levels	1
Maximum output voltage at effect level 1	1300 V _p
Maximum output current	1.25 A
Rated load	300 Ω

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Control characteristics	HPCS
Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.3 to 1.8
Range of power levels	1 to 300 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	460 V _p
Maximum output voltage at effect level 2	525 V _p
Maximum output voltage at effect level 3	700 V _p
Maximum output current	1.5 A
Rated load	500 Ω

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Description of the monopolar coagulation modes:**Fulgurate**

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Waveform	Continuous sinusoidal, modulated
Crest factor	5 to 6.5
Range of power levels	1 to 120 W
Effect levels	1
Maximum output voltage at effect level 1	3600 V _p
Maximum output current	1.5 A
Rated load	500 Ω

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Waveform	Continuous sinusoidal, modulated
Crest factor	5 to 6.5
Range of power levels	1 to 120 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	1270 V _p
Maximum output voltage at effect level 2	1495 V _p
Maximum output voltage at effect level 3	2000 V _p
Maximum output current	1.2 A
Rated load	500 Ω

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Waveform	Continuous sinusoidal, modulated
Crest factor	3.5 to 6
Range of power levels	1 to 120 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	1224 V _p
Maximum output voltage at effect level 2	1519 V _p
Maximum output voltage at effect level 3	2000 V _p
Maximum output current	1.2 A
Rated load	500 Ω

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Waveform	Continuous sinusoidal, unmodulated
Crest factor	1.3 to 1.8
Range of power levels	1 to 200 W
Effect levels	1 to 5
Maximum output voltage at effect level 1	150 V _p
Maximum output voltage at effect level 2	176 V _p
Maximum output voltage at effect level 3	198 V _p
Maximum output voltage at effect level 4	221 V _p
Maximum output voltage at effect level 5	235 V _p
Maximum output current	3.2 A
Rated load	50 Ω

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Waveform	Continuous sinusoidal, modulated
Crest factor	1.3 to 1.8
Range of power levels	1 to 200 W
Effect levels	1
Maximum output voltage at effect level 1	4300 V _p
Maximum output current	0.8 A
Rated load	500 Ω

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Waveform	Continuous sinusoidal, modulated
Crest factor	7.5 to 9
Range of power levels	1 to 120 W
Effect levels	1 to 3
Maximum output voltage at effect level 1	2520 V _p
Maximum output voltage at effect level 2	4300 V _p
Maximum output voltage at effect level 3	4300 V _p
Maximum output current	0.8 A
Rated load	500 Ω

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Recommended compatible neutral electrodes

Olympus Catalog number	Description
WA99512C	Neutral electrode, for newborns, with cable Patient weight < 5 kg
WA99513A	Neutral electrode, for children, without cable Patient weight 5.. 15 kg
WA99513C	Neutral electrode, for children, with cable Patient weight 5.15 kg
WA99514A	Neutral electrode, for adults, without cable Patient weight > 15 kg
WA99514C	Neutral electrode, for adults, with cable Patient weight > 15 kg
WA99516A	Neutral electrode, universal size, without cable Patient weight > 5 kg (adults and children)

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WA99516C	Neutral electrode, universal size, with cable Patient weight > 5 kg (adults and children)
----------	---

Other compatible neutral electrodes

Manufacturer Model	Neutral electrode description
	Patient weight 5...15 kg
3M 8180F	
3M/Medline 9160, 9165	Patient weight 5...15 kg
Bowa 816-141, 818-141, 816-161, 818-161	Patient weight > 15 kg
	Patient weight 5...15 kg
Bowa 816-071, 818-071	
Bowa 816-072, 818-072, 816-092, 818-092, 816-112, 818-112	Patient weight > 5 kg
Conmed 51-7310, 51-7410	Patient weight > 15 kg
	Patient weight 5...15 kg
Conmed 425-2508, 425-2801, 51-7710	
Conmed 410-2000, 410-2100, 410-2200, 410-2400	Patient weight > 2.2 kg
Covidien E7507, E7507-DB, E7509, E7509-B	Patient weight > 13.6 kg
Covidien E7510-25, E7510-25DB	Patient weight 2.7...13.6 kg

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© Copyright 2019 Olympus Winter & Ibe GmbH
© Copyright 2020 Olympus Winter & Ibe GmbH