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
This document is the User manual for iCT motor(WL-1) of Dentium Co., Ltd.

The purpose of notarization is to submit User manual to The Federal Service on Surveillance in Healthcare of Russia.

Dentium Co., Ltd. has the responsibility of any problem which may be caused by this letter of statement.

Date : 23. March. 2023

Company Name : Dentium Co., Ltd

Seung Woo Suh

CEO

6<sup>th</sup> of September, 2021

Approved by

(signature)


Suh Seung Woo

Director

Dentium Co.,Ltd.(ICT Branch)

**INSTRUCTION FOR USE  
FOR THE MEDICAL DEVICE**

**Physiodispenser ICT motor, version WL-1**

**Version 2**

## 1. Medical device name

Physiodispenser ICT motor, version WL-1.

**WL-1 version, consisting of:**

1. Main body – 1 pc.
2. Two-button foot controller – 1 pc.
3. Micromotor – 1 pc.
4. AC power cable – 1 pc.
5. Sterilization plug – 2 pcs.
6. Handpiece rack – 1 pc.
7. Irrigation stanchion – 1 pc.
8. Irrigation tube holder – 15 pcs.
9. Instructions for use – 1 pc.
10. Main body charger – 1 pc.
11. AA battery – 2 pcs.

*Hereinafter referred to as physiodispenser, ICT motor WL-1, WL-1, medical device, device, apparatus.*

## 2. Information about the manufacturer, developer, and authorized representative of the manufacturer

### **Manufacturer:**

Dentium Co., Ltd (ICT Branch)

76, Changnyong-daero 256beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do 16229, (KOREA)

Phone: + 82-70-7098-4027

Fax: + 82-31-281-5561

[www.dentium.co.kr](http://www.dentium.co.kr)

### **Developer**

Dentium Co., Ltd (ICT Branch)

76, Changnyong-daero 256beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do 16229, (KOREA)

Phone: + 82-70-7098-4027

Fax: + 82-31-281-5561

[www.dentium.co.kr](http://www.dentium.co.kr)

## Authorized representative of the manufacturer in the Russian Federation

Limited Liability Company IMPLANT.RU (IMPLANT.RU, LLC)

Bolshoy Gnezdnikovskii alley, 3, Moscow, 125009

Phone/Fax: +7(495)118-76-20, +7(495)941-02-61

## 3. Purpose of the device and potential consumers

### **Purpose of the device**

ICT motor is a device that provides fastening and drive (movement) of handpieces, supply of fluids to the operating area for dental, surgical and implantation procedures.

### **Potential consumer**

Physiodispenser is intended for use by qualified professionals: dentists, dental surgeons, maxillofacial surgeons, implant surgeons.

## 4. Indications, contraindications and possible side effects

Indications for use are dental, surgical (for example, cutting and removing bones, teeth and hard tissue) and implantation procedures (for example, placing an implant or abutment) at the discretion of a

INSTRUCTION FOR USE FOR THE MEDICAL DEVICE

Physiodispenser ICT motor, version WL-1 (version 2)

Page 2 of 22

specialist.

Contraindication is inappropriate use, use with unapproved parts or inappropriate parts.\*

\*The medical device should be used in combination only with handpieces and irrigating tubes, registered in an appropriate manner in the Russian Federation and having the specification listed in the clause 7.1 (table 1 (for handpieces) and table 2 (for irrigating tubes)).

Side effects: ICT motor physiodispenser is an auxiliary device for the operation, therefore does not lead by itself to the development of side effects.

## 5. Precautions

- The ICT motor should be operated by dentists or authorized specialists.
- Read this Instruction for use thoroughly before operation for the first time, in order to avoid unintended incidents or accident.
- If repair or inspection is required, please contact the authorized representative of the manufacturer in the Russian Federation.
- Sterilize all subcomponents that are subject to sterilization before each operation for disinfection. Detailed sterilization instruction is included in this Instruction for use.
- Put the handpiece on the handpiece rack for idling during usage.
- If any abnormal condition occurs, such as unusual noise or overheated components, please stop using the device and contact the authorized representative of the manufacturer in the Russian Federation for support.
- The foregoing warranty shall not apply to defects resulting from unauthorized disassembling or using ingenuine components or parts. In the event of modifications by third parties, the warranty shall become null and void. We recommend using only genuine spare parts for the operation and repair of the ICT motor.
- Do not place similar medical electrical devices near the ICT motor as these devices may disrupt the performance of this device. To mitigate disruptive electromagnetic interference, position the ICT motor away from RF transmitters and other sources of electromagnetic energy.
- Dropping or exerting excessive forces onto the device may cause the device to malfunction or fail to meet the desired performance standard.
- While the burr (cutting tool) is rotating, it should be handled with extra care to avoid injury of the user.
- Do not assemble or disassemble the handpiece with micromotor during operation. Before each usage, try to pull out the handpiece from the micromotor to make sure they are well assembled.
- To avoid risk of electric shock, this equipment must only be connected to supply mains with protective earth.
- Do not modify this equipment without authorization of the manufacturer.
- Do not put your hand in to the irrigation pump gap.

## 6. Equipment classification

- ❖ Electric shock protection and applied part:
  - Class II with internal power supply, type B (control box with control box charger), with internal power supply, type B (two-button foot controller).
- ❖ Method of sterilization and disinfection recommended by the manufacturer:
  - See sub-item Sterilization
- ❖ Waterproof rate of the two-button foot controller: IPX 1
- ❖ Mode of operation
  - Short operating mode

## 7. Description and design of the device

### 7.1. Device description

The medical device «Physiodispenser ICT motor» is a device, the essence of which is to fix and insert handpieces, deliver solutions to the operation area for cutting bones, placing an implant or abutment, and other functions (see indications for use). The design of the motor allows for variable speed (rpm) and torque, and a selection of available options and/or cooling water supply based on the surgical phase. Various ratios of the final value of the speed of rotation and the torque are provided in the engine in accordance with the gear ratio during mechanical reduction.


The device is used in conjunction with dental handpieces and disposable sterile irrigation tubes. The handpiece used with the micromotor should be registered in the territory of the country where the device will be used in the established order and have the characteristics presented in table 1.

**Table 1. Parameters of handpieces to be used with physiodispenser.**

| Gear ratio | Maximum torque |
|------------|----------------|
| 20:1       | 5-70 N-cm      |


Sterile disposable irrigation tubes must also be registered in the country where they will be used according to the established procedure. The dimensional characteristics of the irrigation tubes compatible with the device are shown in Table 2.

**Figure 1. Sterile disposable irrigation tubes.**

**Table 2. Irrigation tube dimensions**

| Segment | Size   |
|---------|--------|
| A       | 123.5  |
| B       | 240    |
| C       | 93 ± 5 |
| D       | 2.000  |
| (1)     | Ø5.5   |
| (2)     | Ø2.6   |

Physiodispenser is intended for repeated use and is a device with a short operating mode. The physiodispenser has no contact with the patient and the operator. The doctor works with gloves.


### 7.2. Package contents

#### ① Main body


**Figure 2. Main body. General view.**

Front side.

Back side.

**Table 3. Main body description.**

INSTRUCTION FOR USE FOR THE MEDICAL DEVICE  
Physiodispenser ICT motor, version WL-1 (version 2)

| No. | Name                               | Description                                                                                                                                         |
|-----|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| (a) | LCD Panel                          | Displays setting value and operation status                                                                                                         |
| (b) | Mode Change Arrow (Left)           | The button to move the user program ("STEP") up and down (Countersinking, Implant Insertion or Tapping, Cover & Healing connection, Removal, Spare) |
| (c) | Mode Change Arrow (Right)          |                                                                                                                                                     |
| (d) | Micromotor Connector               | Micromotor connection connector                                                                                                                     |
| (e) | Irrigation Pump Module             | Irrigation solution pumping module with immersion tube mounting                                                                                     |
| (f) | Irrigation Module Open/Close Lever | Lever used when mounting irrigation tube to immersion module                                                                                        |
| (g) | Set Value Up/Down                  | Switch to change setting value after selecting mode (Up / Down)                                                                                     |
| (h) | Irrigation pump stanchion hole     | The hole to install an irrigation stanchion                                                                                                         |
| (i) | Power Switch                       | Power ON/OFF switch                                                                                                                                 |
| (j) | DC Jack                            | Connect charger adapter when using equipment or charging DC power jack                                                                              |

## LCD Panel

Figure 3. LCD Panel.


Table 4. Description of the LCD Panel

| No.                                                                                 | Name                                                                                  | Description                                                                                                                                                                                                                                                 |         |         |                                                                                     |                                                                                       |       |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------|---------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|-------|
| 1                                                                                   | Show remaining battery                                                                | 4 stages                                                                                                                                                                |         |         |                                                                                     |                                                                                       |       |
| 2                                                                                   | Change and show the program                                                           | 0~7 stages                                                                                                                                                             |         |         |                                                                                     |                                                                                       |       |
| 3                                                                                   | Set and show micromotor rotation speed (rpm)                                          | 10~2,000 rpm                                                                                                                                                                                                                                                |         |         |                                                                                     |                                                                                       |       |
| 4                                                                                   | Set and show micromotor torque                                                        | 5~70 Ncm                                                                                                                                                                                                                                                    |         |         |                                                                                     |                                                                                       |       |
| 5                                                                                   | Show remaining battery of foot controller                                             | 4 stages                                                                                                                                                                |         |         |                                                                                     |                                                                                       |       |
| 6                                                                                   | Show wireless communication                                                           |                                                                                                                                                                        |         |         |                                                                                     |                                                                                       |       |
| 7                                                                                   | Change and show mode                                                                  | 0~7                                                                                                                                                                    |         |         |                                                                                     |                                                                                       |       |
| 8                                                                                   | Set and show forward and reverse operation of micromotor                              | F/R                                                                                                                                                                                                                                                         |         |         |                                                                                     |                                                                                       |       |
|                                                                                     |                                                                                       | <table><tr><th>Forward</th><th>Reverse</th></tr><tr><td></td><td></td></tr></table> | Forward | Reverse |  |  |       |
| Forward                                                                             | Reverse                                                                               |                                                                                                                                                                                                                                                             |         |         |                                                                                     |                                                                                       |       |
|  |  |                                                                                                                                                                                                                                                             |         |         |                                                                                     |                                                                                       |       |
| 9                                                                                   | Change and show irrigation (water supply)                                             | 0~3 ((unit: mL/min)) <table><tr><td>Water</td><td>OFF</td><td>40 mL</td><td>50 mL</td><td>60 mL</td></tr></table>                                                                                                                                           | Water   | OFF     | 40 mL                                                                               | 50 mL                                                                                 | 60 mL |
| Water                                                                               | OFF                                                                                   | 40 mL                                                                                                                                                                                                                                                       | 50 mL   | 60 mL   |                                                                                     |                                                                                       |       |

|    |                                               |                  |                                                                                   |                                                                                     |                                                                                     |                                                                                     |
|----|-----------------------------------------------|------------------|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|    |                                               | Display          |  |  |  |  |
| 10 | Adjust and show brightness of handpiece optic |                  | 0~3                                                                               |                                                                                     |                                                                                     |                                                                                     |
|    |                                               | Bright-ness      | OFF                                                                               | Step 1                                                                              | Step 2                                                                              | Step 3                                                                              |
|    |                                               | Display          |  |  |  |  |
| 11 | Change and show gear ratio                    | 20:1, 27:1, 32:1 |                                                                                   |                                                                                     |                                                                                     |                                                                                     |
| 12 | Show handpiece calibration                    | Calibration      |                                                                                   |                                                                                     |                                                                                     |                                                                                     |

## ② Two-button foot controller

Figure 4. Two-button foot controller (general view and side view).


Figure 5. Two-button foot controller (bottom view).


Table 5. Description of the two-button foot controller.

| No. | Name                       | Description                                                                                      |
|-----|----------------------------|--------------------------------------------------------------------------------------------------|
| ①   | Two-button foot controller | Micro-motor operation by pressing foot switch / Change speed                                     |
| ②   | Mode Change Switch         | When the button is pressed, the program is operated step by step (0-1-2-3-4-5-6-7) / Change mode |
| ③   | Power Switch               | Power On/Off switch for operation                                                                |

|     |                           |                                                  |  |
|-----|---------------------------|--------------------------------------------------|--|
|     |                           | OFF/ON                                           |  |
| (d) | Battery Protection Cover* | Battery mount for foot controller operation      |  |
| (e) | Battery Holder            | The area where the battery is attached and fixed |  |

\* the set contains two AA batteries

### ③ Micromotor

Figure 6. Micromotor (general view).



Table 6. Description of Micromotor

| No.    | Name                           | Description                                                                                                      |
|--------|--------------------------------|------------------------------------------------------------------------------------------------------------------|
| (a)    | Micromotor Adapter             | The part connected with the handpiece.                                                                           |
| (b)    | Micromotor                     | A BLDC motor to give torque to contra angle handpiece. A rpm and torque are controlled by signal from main body. |
| (c)(d) | Micromotor cable and connector | A part to connect with micromotor in main body.                                                                  |

### ④ AC power cable, sterilization plug, handpiece rack, irrigation stanchion, irrigation tube holder, adapter

Table 7. AC power cable, sterilization plug, handpiece rack, irrigation stanchion, irrigation tube holder description

| Name           | Photo                                                                                | Description                                                                                         |
|----------------|--------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------|
| AC power cable |  | Designed to transmit electricity to the main body and is connected to the socket on its rear panel. |

|                        |                                                                                     |                                                                                                                  |
|------------------------|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------|
| Sterilization plug     |    | Fits onto the handpiece adapter during sterilization. Designed to protect the micromotor from steam penetration. |
| Handpiece rack         |    | Stand for positioning the handpiece of the micromotor between use.                                               |
| Irrigation Stanchion   |   | Designed to fix the package with the irrigation solution.                                                        |
| Irrigation tube holder |  | Designed to hold the irrigation tubes, installed on the micromotor cable and on the irrigation tube.             |

## ⑤ Main body charger.

Figure 7. Main body charger (general view)


Table 8. Description of the main body charger

| Name    | Description                                                       |
|---------|-------------------------------------------------------------------|
| Adapter | The Adapter is designed to supply power to the main body (42 VDC) |
| Cord    | Designed to connect the Adapter to the main body                  |

## 8. Operation and application procedure

### 8.1. Operation

The design of the micromotor allows to adjust the rotational speed (rpm) and torque of the dental handpiece, to select the available options and/or the supply of cooling water, taking into account the stage of the surgical procedure. Different ratios of the final value of the rotational speed and the torque in the motor are provided in accordance with the gear ratio during mechanical reduction (reduction of the handpiece revolution number).

- When the system power switch is on, it converts 230 VAC to 36 VDC.
- The converted 36V power supply is fed to the system in a brushless DC motor drive, separated from the circuit power (5VDC) and the motor power.
- When the power is on, the system is activated, the system checks and LCD screen is displayed.
- When the foot switch is pressed during the procedure, the program sends a control signal to the micromotor according to the analog signal of the foot switch, then the micromotor starts rotating. At this time, the irrigation pump also starts to operate at the set speed.
- Main control: LCD screen + footswitch + motor signal control.
- Brushless DC motor drive: it drives the micromotor and the irrigation motor

## 8.2. Application procedure

### 8.2.2. Product installation

#### Power Connection

- Plug the power connector into the socket on the back of the main body to supply power.

- Use only the charger provided by manufacturer. If any charge other than provided is used, it may cause the fire or electric shock.
- Dry your hands completely before connecting the power cable.
- This product has the free voltage function (AC 100-240 V~, 50-60 Hz).
- To avoid risk of electric shock, this equipment must only be connected to supply mains with protective earth.

### 8.2.3. Irrigation system installation

- Place the irrigation solution stand-holder in the hole in the control unit box.
- Install the irrigation tube:
  - Open the cover to install the tube by moving the OPEN/CLOSE lever to the OPEN position.
  - Insert the irrigation tube into the adapter with both hands. Place the irrigation tube exactly in the middle of the roller to avoid kinking or crushing. Once the irrigation tube is properly seated in the tube guide, move the OPEN/CLOSE lever to the CLOSE position.
- Attach the handpiece to the irrigation tube. Refer to the handpiece supplier's user manual for detailed instructions.

### 8.2.4. Connecting the components

#### Connecting and disconnecting of micromotor

- Place the [↔] mark on the micromotor cord so that it is on top and insert the micromotor plug into the micromotor socket on the front panel of the control unit. A click indicates correct connection. Take the micromotor plug and pull towards you to disconnect the micromotor. It is imperative to turn off the power supply to the main unit when connecting and disconnecting the plugs.


The LCD screen will show an error notification if the micromotor is disconnected or incorrectly connected. The connection should be checked for normal use.

### Connection of two-button foot controller

- As the wireless module is applied to ICT motor foot switch, if the power of the foot Controller and mail unit is on, wireless communication is connected automatically.
- Turn the power on by pressing the power switch toward  direction.
- Two-button foot controller serves to operate or stop the equipment. To stop the operation, take the foot off two-button foot controller.

### 8.2.3. LCD screen application procedure

#### Turning Power Switch On

Turn on the power by pressing power switch located rear of the main unit toward  direction.

**Table 9. Control panel functions**

| Button                                                                              | Detailed Function                                                                                                                                                                                                                                                                                                                                                                               |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | When pressing this button, green lamp at the bottom of setting mode in LCD display is turned on and move by one stage whenever it is pressed (Program → F/R → Irrigation → Lamp → Gear Ratio → Calibration → Torque → RPM)                                                                                                                                                                      |
|  | When pressing this button, green lamp at the bottom of setting mode in LCD display is turned on and move by one stage whenever it is pressed (Program → RPM → Torque → F/R → Irrigation → Lamp → Gear Ratio → Calibration)                                                                                                                                                                      |
|  | <ul style="list-style-type: none"> <li>- When rotating jog dial after selecting the mode, the value of relevant function can be changed. It is used by setting normal direction/reverse direction.</li> <li>- When you want to save the settings in memory after completing the setting by each mode, if the button is pressed for 3 seconds, the setting is saved in relevant mode.</li> </ul> |

## Program Selection and Change

ICT motor provides you with the optimized program solution for each STEP of protocol for the Dentium implant installation as default from the factory. This program can be changed by the user according to preferences.


**Table 10. Program (Step) Information\***

| STEP | Ratio | Process                              | rpm  | Torque | Pump | Direction |
|------|-------|--------------------------------------|------|--------|------|-----------|
| 0    | 20:1  | Drilling                             | 0    | 5      |      |           |
| 1    | 20:1  | Drilling                             | 1000 | 35     |      |           |
| 2    | 20:1  | Countersinking                       | 300  | 35     |      |           |
| 3    | 20:1  | Implant Insertion or Tapping         | 20   | 35     |      |           |
| 4    | 20:1  | Cover screw & Healing cap connection | 100  | 5      |      |           |
| 5    | 20:1  | Removal                              | 20   | 55     |      |           |
| 6    | 20:1  | Spare                                | 20   | 5      |      |           |
| 7    | 20:1  | Spare                                | 0    | 0      |      |           |

\*For your reference, above are not the standard guidance for all implant system.

## Changing a STEP

- When moving the program, if pressing switch, the blue selection lamp is turned on and if pressing the switch, move to left and right by step. And when the blue selection lamp is located below the Program, set desired stage by rotating jog dial.
- For the method to store the value for the program, if pressing jog dial for 3 seconds after completing the setting, the value will be saved in set program together with beep sound.

### Changing rotation speed (rpm)

- When moving the program, if pressing   switch, the blue selection lamp is turned on and if pressing the switch, move to left and right by step. And when the blue selection lamp is located below RPM, set desired stage by rotating jog dial.
- For the method to store the value for the program, if pressing jog dial for 3 seconds after completing the setting, the value will be saved in set program together with beep sound.

- The unit to be changed according to rotational speed are shown in below.

**Table 11. Increment value steps.**

| Setting rpm | Increment rpm |
|-------------|---------------|
| 15-50       | 5             |
| 60-100      | 10            |
| 200-2000    | 100           |

- The performance may be unstable slightly ( $\pm 10\%$ ) under the low speed of 25 RPM or less.

### Changing Torque

- When moving the program, if pressing   switch, the blue selection lamp is turned on and if pressing the switch, move to left and right by step. And when the blue selection lamp is located below N·cm, set desired stage by rotating jog dial.
- For the method to store the value for the program, if pressing jog dial for 3 seconds after completing the setting, the value will be saved in set program together with beep sound.

- Torque value is increased or decreased by 5 N·cm unit. The performance may be unstable slightly ( $\pm 10\%$ ) under the torque of 55 N·cm.

### Changing Direction

- When moving the program, if pressing   switch, the blue selection lamp is turned on and if pressing the switch, move to left and right by step. And when the blue selection lamp is located below , set desired stage by rotating jog dial.


- For the method to store the value for the program, if pressing jog dial for 3 seconds after completing the setting, the value will be saved in set program together with beep sound.

- When setting to counter-clock wise ("R"), long alarm sound is made continuously and after operating the motor, it is converted to short alarm sound and lasts continuously.
- All the set values can be changed and manipulated only when the motor operation is completely stopped.
- "F" means "Forward" and "R" means "Reverse".

### Changing Gear Ratio

- When moving the program, if pressing   switch, the blue selection lamp is turned on and if pressing the switch, move to left and right by step. And when the blue selection lamp is located below 20:1, set desired stage by rotating jog dial.
- For the method to store the value for the program, if pressing jog dial for 3 seconds after completing the setting, the value will be saved in set program together with beep sound.

See the table 12 for the speed and torque range to be used with the change in gear ratio.

**Table 12. Value ranges.**

| rpm  | rpm     | Torque    |
|------|---------|-----------|
| 20:1 | 10-2000 | 5-70 N·cm |

### Changing Irrigation Volume

- When moving the program, if pressing   switch, the blue selection lamp is turned on and if pressing the switch, move to left and right by step. And when the blue selection lamp is

located below , set desired stage by rotating jog dial.

- For the method to store the value for the program, if pressing jog dial for 3 seconds after completing the setting, the value will be saved in set program together with beep sound. Irrigation can be changed with 4 stages and for the output by each stage, see below.

**Table 13. Irrigation volume ranges.**

| Stage                                                                             | Irrigation volume (Max.) mL/min |
|-----------------------------------------------------------------------------------|---------------------------------|
|  | OFF                             |
|  | 40 mL                           |
|  | 50 mL                           |
|  | 60 mL                           |

- The device makes short beep sounds when the actual torque value reaches 90 % of the target torque value. Once torque value reaches 100 % of the target torque, the device makes long beep sound and the Micromotor will stop in 0.5 second.
- The device will make short beep sounds if the handpiece is working in reverse direction.
- The parameters can be adjustable (STEP change) only if the device has come to complete stop.

## 9. Technical Specification

**Table 14. Main technical characteristics**

| Parameters                                        | Specification                                                                                                                                           |
|---------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| Input power voltage and frequency                 | 100-240 V ~ 50/60 Hz, 130-160 V·A                                                                                                                       |
| Electric shock protection and applied part        | Class II with internal power supply, type B (control box with control box charger)<br>With internal power supply, type B (two-button foot controller)   |
| Handpiece speed and handpiece torque values range | 20-2,000 rpm                                                                                                                                            |
| Handpiece torque range                            | 5-70 N·cm<br>(based on gear ratio 20:1)                                                                                                                 |
| Main body input rated power                       | 42 VDC, 1.5 A (when connected to the charger of the main unit)<br>36.5 VDC, 2,750 mAh (lithium-ion battery only)                                        |
| Adapter input rated power                         | 100-240 V ~ 50/60 Hz, max. 1.5 A                                                                                                                        |
| Adapter output rated power                        | 42 VDC, 1.5 A                                                                                                                                           |
| Two-button foot controller input rated power      | 3 VDC                                                                                                                                                   |
| Wireless communication specifications             | Frequency range: 2,405-2,480 MHz<br>Number of channels: 16<br>Channel Separation: 5 MHz<br>RF output power: 8.6 dBm Equivalent isotropic Radiated Power |
| Waterproof rate of the two-button foot controller | IPX 1                                                                                                                                                   |
| Irrigation speed                                  | Level 1: 40 ml/min<br>Level 2: 50 ml/min<br>Level 3: 60 ml/min                                                                                          |
| LCD screen                                        | Longitudinal profile LCD                                                                                                                                |
| Mode of operation                                 | Short operating mode.<br>Maximum activation time: 3 min<br>Minimum deactivation time: 10 min                                                            |
| Sound pressure level                              | Not more than 80 dBA                                                                                                                                    |

### 9.1. Information about electromagnetic compatibility

- This equipment is intended for use in the electromagnetic environment specified below.
- The customer or the user of this equipment should assure it is used in such an environment.

**Table 15. Electromagnetic compatibility**

| Emissions test                         | Compliance               | Electromagnetic environment-guidance                                                                                                                                                                                                                                                                                                                                                                                    |
|----------------------------------------|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Radiated disturbance                   | CISPR 11 Group 1 Class A | The EMISSIONS characteristics of ICT motor make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment. |
| Mains terminal disturbance voltage     | CISPR 11 Group 1 Class A | The ICT motor is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded:                                                                                                                       |
| Harmonic emissions                     | Class A                  |                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Voltage fluctuations flicker emissions | Complies                 | Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the ICT motor or shielding the location.                                                                                                  |

### 9.2. Electromagnetic immunity

- This equipment is intended for use in the electromagnetic environment specified below.
- The customer or the user of this equipment should assure it is used in such an environment.

**Table 16. Electromagnetic immunity**

| Immunity test                                                  | IEC 60601 test level                                           | Conformance level                                              | Electromagnetic environment-guidance                                                                                                                        |
|----------------------------------------------------------------|----------------------------------------------------------------|----------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge (ESD) IEC 61000-4-2                    | Discharge by 8 kV direct contact<br>15 kV of Air-gap discharge | Discharge by 8 kV direct contact<br>15 kV of Air-gap discharge | Floors should be made of wood or concrete or have ceramic tiles. When the floor is made of synthetic material, the relative humidity must be at least 30 %. |
| Radiated RF<br>Electromagnetic Field immunity<br>IEC 61000-4-3 | 3 V/m<br>80 MHz-2.7 GHz<br>80 % AM at 1 kHz                    | 3 V/m<br>80 MHz-2.7 GHz<br>80 % AM at 1 kHz                    | The ICT motor is suitable to use in professional healthcare environment.                                                                                    |


|                                                                                                     |                                                                                                                                                                 |                                                                                                                                                                 |                                                                                                                                                                                                                                                                                                                                |
|-----------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunity to Proximity Fields from RF wireless Communications Equipment<br>Table 9 in IEC 60601-1-2  | 28 V/m Max.<br>385-5785 MHz in according to table 9                                                                                                             | 28 V/m Max.<br>385-5785 MHz in according to table 9                                                                                                             | RF communication equipment is used on closer than 30 cm to any part of the ICT motor, including cables specified by Dentium.                                                                                                                                                                                                   |
| Electrical fast transient/burst IEC 61000-4-4                                                       | $\pm 2$ kV for power supply lines<br>$\pm 1$ kV for I/O lines<br>100 kHz repetition frequency                                                                   | $\pm 2$ kV for power supply lines<br>$\pm 1$ kV for I/O lines<br>100 kHz repetition frequency                                                                   | The quality of supplied power should be suitable for general business site or hospital environment.                                                                                                                                                                                                                            |
| Surge IEC 61000-4-5                                                                                 | $\pm 1$ kV line to line<br>$\pm 2$ kV line to earth                                                                                                             | $\pm 1$ kV line to line<br>$\pm 2$ kV line to earth                                                                                                             | The quality of supplied power should be suitable for general business site or hospital environment.                                                                                                                                                                                                                            |
| Immunity to conducted disturbances induced by RF fields IEC 61000-4-6                               | 3 V<br>0.15-80 MHz<br>6 V in ISM bands between 0.15 and 80 MHz<br>80 % AM at 1 kHz<br><br>Power supply line & I/O lines                                         | 3 V<br>0.15-80 MHz<br>6 V in ISM bands between 0.15 and 80 MHz<br>80 % AM at 1 kHz<br><br>Power supply line & I/O lines                                         | The strength of RF field in the frequency range higher than 150 kHz~80 MHz, the strength of the RF field is smaller than 3 V.                                                                                                                                                                                                  |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | 0 % UT: 0.5 cycle<br>At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°<br>0 % UT; 1 cycle and 70 % UT; 25/30 Cycles Single phase: at 0°<br>0 % UT; 250/300 cycle | 0 % UT: 0.5 cycle<br>At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°<br>0 % UT; 1 cycle and 70 % UT; 25/30 Cycles Single phase: at 0°<br>0 % UT; 250/300 cycle | The quality of supplied power should be suitable for general business site or hospital environment.<br>For the user to operate the equipment continuously even the electric power supply is interrupted, it is recommended that the uninterruptable power supply device (UPS) or battery is prepared.                          |
| Magnetic field of supply frequency (50/60 Hz) IEC 61000-4-8                                         | 30 A/m                                                                                                                                                          | 30 A/m                                                                                                                                                          | If image distortion occurs, it may be necessary to position the Model ICT motor image intensifier further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low. |

### 9.3. Protection Distance

- Portable and mobile HF communication equipment should be operated at recommended protection distance from this equipment.
- This equipment is manufactured to operate in the electromagnetic environment where the emitted HF disturbance is controlled.
  - RF communication equipment is used not closer than 30 cm to any part of the ICT motor including cables specified by Dentium Co.Ltd. (ICT Branch).

## 10. Manufacturer's package and labelling

Table 17. Symbols and designations

| Symbol                                                                              | Description                                                                                |
|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|
|    | Type B applied part                                                                        |
|    | Refer to the instructions for use for important information, namely warnings and cautions. |
|    | Caution! Moving machinery                                                                  |
|    | Fuse                                                                                       |
|    | Foot controller (two-button foot controller)                                               |
|    | Do not dispose of with household waste                                                     |
|  | Refer to operation manual/leaflet                                                          |
|  | Caution!                                                                                   |
|  | Note                                                                                       |
|  | Compliance with European standards                                                         |
|  | Medical device manufacturer                                                                |
|  | Date of manufacture of the medical device                                                  |
|  | Authorized Representative in the EU                                                        |
|  | Catalog number                                                                             |
|  | Serial number                                                                              |
|  | Water protection degree                                                                    |
|  | Autoclaving (up to a maximum temperature of 135 °C)                                        |

|                                                                                     |                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------|
|    | Autoclaving (up to a maximum temperature of 134 °C)   |
|    | Top (correct vertical position of the load)           |
|    | Inadmissibility of shock loading                      |
|    | Limit on the number of tiers in a stack               |
|    | Moisture control                                      |
|    | "On" position (state) of the device part power supply |
|  | "On" position (state) of the device part power supply |
|  | Power plug                                            |
|  | Alternating current                                   |
|  | Direct current                                        |
|  | Class II product                                      |

## 11. Commissioning, maintenance and repair

### Commissioning conditions

- Cables should be checked for peeling of the insulating cover and stored in a roll.
- Stop further use of the device and contact an authorized representative of the manufacturer for service in case of technical problems, such as noise or overheating during operation.
- Dispose of any part of electronic equipment in an environmentally friendly way.
- It is unacceptable to locate the equipment or any part of it in a place exposed to water or moisture.
- Locate the equipment in a location free from exposure to temperature, ventilation, sunlight, dust, salt, etc.
- Avoid placing the equipment on an inclined surface or subjecting it to shock.
- Do not store equipment near chemicals or gases.
- Disconnect power to the equipment and apply to professional assistance if water or liquid enters the equipment.

**Table 18. Checklist of daily check parameters**

| Description                | Check item                                                                                                                                                                                                                                                                                                                          |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Main body                  | <ol style="list-style-type: none"> <li>1. Is the I/O power switch of unused equipment in the "⊙" position.</li> <li>2. Are any switches on the control panel damaged.</li> <li>3. Whether the colors and data are clearly displayed on the screen.</li> </ol>                                                                       |
| Micromotor                 | <ol style="list-style-type: none"> <li>1. Whether the micromotor cables are damaged.</li> <li>2. Whether the connector of the micromotor, which is used to connect to the main body, is not damaged.</li> <li>3. Is the connection of any part of the micromotor assembly loose.</li> </ol>                                         |
| Two-button foot controller | <ol style="list-style-type: none"> <li>1. Whether the two-button foot controller power switch is operating as intended.</li> <li>2. Whether the two-button foot controller moves smoothly when you press it.</li> <li>3. Whether the button for changing the mode of the two-button foot controller is pressed smoothly.</li> </ol> |

The battery used in the main body has a service life of about 2 years. Please contact an authorized representative of the manufacturer if the battery level displayed on the screen is low during a routine check.

### Troubleshooting

- Inspection of the device should be carried out under conditions that exclude the presence of water or moisture.
- Troubleshooting should be done after disconnecting the power supply.
- A list of simple malfunctions that can be corrected by the user is given below. You should contact the authorized representative of the manufacturer in the Russian Federation if none of the suggested options helped to fix the problem.

**Table 19. Possible reasons of malfunctions and solutions.**

| Malfunction                  | Reason and solution |                                                                                                                                                                                                                                                                                                                                                                                                 |
|------------------------------|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| LCD does not turn on         | Reason              | Lack of power supply                                                                                                                                                                                                                                                                                                                                                                            |
|                              | Solution            | <ol style="list-style-type: none"> <li>1. Unplug the main body's power connector from the power outlet and reconnect it.</li> <li>2. Turn the main power switch to the "⊙" position.</li> <li>3. Use a battery charged with a charger.</li> <li>4. If the performance of the battery deteriorates, connect the charger.</li> </ol>                                                              |
| Control button does not work | Reason              | <ol style="list-style-type: none"> <li>1. The button of the running equipment does not work.</li> <li>2. The control button is defective.</li> </ol>                                                                                                                                                                                                                                            |
|                              | Solution            | <ol style="list-style-type: none"> <li>1. Wait until the operation is complete and press the button.</li> <li>2. Contact your nearest Dentium branch for service.</li> </ol>                                                                                                                                                                                                                    |
| Handpiece does not work      | Reason              | <ol style="list-style-type: none"> <li>1. The handpiece does not work because the power supply of the two-button foot controller is turned off.</li> <li>2. The handpiece does not work because it is dirty.</li> <li>3. The tip does not work due to wear or breakage of the internal gear.</li> <li>4. The handpiece does not work, although the two-button foot controller is on.</li> </ol> |

|                                                      |          |                                                                                                                                                                                                                                                                                                                                                                                                   |
|------------------------------------------------------|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                      | Solution | <ol style="list-style-type: none"> <li>1. Check if the handpiece is working by turning on the power supply to the two-button foot controller.</li> <li>2. Read the operation manual for the handpiece.</li> <li>3. Connect another handpiece to the equipment and check if it works.</li> <li>4. Contact the authorized representative of the manufacturer in the Russian Federation..</li> </ol> |
| Irrigation pump does not work or is not working well | Reason   | <ol style="list-style-type: none"> <li>1. A tube is damaged or a tip is blocked up.</li> <li>2. The irrigation tube stopper is blocked.</li> <li>3. Weak irrigation mode is set.</li> </ol>                                                                                                                                                                                                       |
|                                                      | Solution | <ol style="list-style-type: none"> <li>1. Replace a tube and check tip attachment. Handpiece cleaning should be carried out in accordance with the instructions in the manufacturer's instruction manual.</li> <li>2. Unlock the irrigation tube stopper.</li> <li>3. Adjust the value of the irrigation volume on the screen (1-3).</li> </ol>                                                   |

## 12. Transportation, storage and operation conditions

**Storage.** Devices should be stored in well-ventilated areas away from heat, sunlight, direct exposure to light sources and identical electrical equipment.

**Transportation.** Transportation of a device packed in a shipping container should be carried out by all types of covered vehicles in compliance with the precautions indicated on the shipping container and in accordance with the rules for the carriage of goods in force for this type of transport.

**Operation.** Devices are used in medical institutions, in conditions requiring the protection of the hands of a medical worker, protective clothing and protection of the thyroid gland for the patient when using the devices as intended.

**Table 20. Operation conditions, transportation and storage**

|                     |                                                                             |
|---------------------|-----------------------------------------------------------------------------|
| Ambient temperature | 10 °C to 35 °C (operation)<br>От 0 °C до 60 °C (transportation and storage) |
| Relative humidity   | Up to 80 % (operation)<br>Up to 90 % (transportation and storage)           |
| Atmosphere pressure | 800-1060 hPa (operation)<br>700-1060 hPa (transportation and storage)       |

## 13. Cleaning, disinfection and sterilization of the device

### Cleaning and disinfection

Cleaning and disinfection of the micromotor and handpiece rack are necessary after each patient. The specified components should be wiped with a clean cloth soaked in a product containing at least 70 % alcohol.

The rest of the components are recommended to be periodically cleaned with a product containing at least 70 % alcohol.


### Sterilization

High pressure steam sterilization is intended for the micromotor and handpiece rack after cleaning. Working conditions for high pressure steam sterilization: 132 °C, sterilization duration: 15 minutes. Maximum number of sterilization sessions: 100 times.

Close the hole in the handpiece connector with a sterilization plug as shown below to prevent steam from entering the micromotor as shown below.

**Figure 8. Closure of the handpiece connector with a sterilization plug.**

#### **14. Disposal**

Disposal of the device should be carried out in accordance with federal and regional regulations.

Separation into parts and components requiring different disposal methods should be done in preparation for disposal:

- units and parts containing electronic components;
- rechargeable batteries;
- mechanical parts of the apparatus.

The micromotor and handpiece should be disinfected before disposal in accordance with the instructions for use.

It is allowed to dispose of parts of the apparatus (after disinfection) as epidemiologically safe medical waste, close to solid household waste in composition.

Caution: Batteries should be disposed of separately from household waste at a collection point for the recycling of electrical and electronic equipment.

The medical organization carries out only collection, disinfection treatment and temporary storage of the above-mentioned waste, until their subsequent removal.

Transportation, disassembly, sorting of medical equipment, their recycling (disposal) and destruction are carried out by specialized organizations licensed to dispose of medical equipment.

#### **15. Shelf life and warranty period**

The shelf life for this type of device is not applicable.

The warranty period of the medical device is 1 year after installation.

#### **16. List of international standards that the medical device complies with**

EN ISO 13485, EN ISO 14971, EN 60601-1, IEC 80601-2-60, EN 1640, EN 60601-1-2, EN 60601-1-6, EN 62366, EN 62304, ISO 17665-1, ISO 14457, EN ISO 3964, EN ISO 15223-1, EN ISO 9687, EN 1041, EN ISO 14155, MEDDEV 2.12\_1, rev. 8, MEDDEV 2.7.1 rev. 4.

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Hyo Jeong Lee -----  
-----  
attorney-in-fact of  
Seung Woo Suh -----  
CEO -----

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Dentium Co., Ltd -----  
-----  
appeared before me and admitted  
said principal`s subscription to  
the attached -----  
Letter of Statement. -----  
(User manual) -----

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This is hereby attested on this  
23rd day of Mar. 2023 at this office.

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This office has been authorized by the  
Minister of Justice, the Republic of  
Korea, to act as Notary Public Since  
7, Feb. 2020 Under Law No.138.

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This document is the User manual for iCT motor(WH-1) of Dentium Co., Ltd.  
The purpose of notarization is to submit User manual to the Federal Service on  
Surveillance in Healthcare of Russia.

Dentium Co., Ltd. has the responsibility of any problem which may be caused by this  
letter of statement.

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**Physiodispenser ICT motor, version WH-1**

**Version 2**

## 1. Medical device name

Physiodispenser ICT motor, version WH-1.

### WH-1 version, consisting of:

1. Main body – 1 pc.
2. Two-button foot controller – 1 pc.
3. Micromotor – 1 pc.
4. AC power cable – 1 pc.
5. Sterilization plug – 2 pcs.
6. Handpiece rack – 1 pc.
7. Irrigation stanchion – 1 pc.
8. Irrigation tube holder – 15 pcs.
9. Instructions for use – 1 pc.

*Hereinafter referred to as physiodispenser, ICT motor WH-1, medical device, device, apparatus.*

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
Fax: + 82-31-281-5561

[www.dentium.co.kr](http://www.dentium.co.kr)

### Authorized representative of the manufacturer in the Russian Federation

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Indications for use are dental, surgical (for example, cutting and removing bones, teeth and hard tissue) and implantation procedures (for example, placing an implant or abutment) at the discretion of a specialist.

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Physiodispenser ICT motor, version WH-1 (version 2)

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## 6. Equipment classification

- ❖ Electric shock protection and applied part:
  - Class 1, type B
- ❖ Method of sterilization and disinfection recommended by the manufacturer:
  - See sub-item Sterilization
- ❖ Waterproof rate of the two-button foot controller:
  - IPX 1
- ❖ Mode of operation
  - Short operating mode

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The device is used in conjunction with dental handpieces and disposable sterile irrigation tubes.


The handpiece used with the micromotor should be registered in the territory of the country where the device will be used in the established order and have the characteristics presented in table 1.

**Table 1. Parameters of handpieces to be used with physiodispenser**

| Gear ratio | Maximum torque |
|------------|----------------|
| 20:1       | 5-70 N-cm      |

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
**Figure 1. Sterile disposable irrigation tubes**

**Table 2. Irrigation tube dimensions**

| Segment | Size   |
|---------|--------|
| A       | 123.5  |
| B       | 240    |
| C       | 93 ± 5 |
| D       | 2.000  |
| (1)     | Ø5.5   |
| (2)     | Ø2.6   |

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


Back side.

**Table 3. Description of the main body**


| Name                              | Description                                                                                                                                         |
|-----------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Program up/down                   | The button to move the user program ("STEP") up and down (Countersinking, Implant Insertion or Tapping, Cover & Healing connection, Removal, Spare) |
| Data up/down                      | The button to increase/decrease the data (direction, torque, ratio)                                                                                 |
| Mode/Memory                       | The button to activate data to save the changed data                                                                                                |
| Handpiece connector               | The connector to connect handpiece connector                                                                                                        |
| Irrigation pump stanchion hole    | The hole to install an irrigation stanchion                                                                                                         |
| Open/close lever                  | The lever to open/close an irrigation pump cover                                                                                                    |
| Power switch                      | The button to turn on/off                                                                                                                           |
| Power socket                      | The socket to connect AC power cable                                                                                                                |
| Two-button foot controller socket | The socket to connect two-button foot controller connector                                                                                          |

**Figure 3. Inside view of the main body**

## LCD display description

Figure 4. LCD display

- A) Displays the current STEP (STEP 1~6).
- B) Displays the rpm currently set. (It displays the real-time rpm while the motor runs.)
- C) Displays the torque currently set. (It displays the real-time torque while the motor runs.)
- D) Displays the rpm gauge setting as the ratio graphic.
- E) Displays the rotational direction of the handpiece. (The letters F and R stand for "FORWARD" and "REVERSE" respectively.)
- F) Displays the gear reduction ratio of the handpiece.
- G) Displays the irrigation flow rate.

## ② Two-button foot controller

Figure 5. Two-button foot controller (general view)



Table 4. Description of the two-button foot controller

| Name                              | Description                                                                                |
|-----------------------------------|--------------------------------------------------------------------------------------------|
| Foot switch                       | A switch to operate micromotor. The rpm is changed by foot switch depth.                   |
| Program (STEP) shift              | A button to move program (STEP 1~6).                                                       |
| Two-button foot controller socket | A part to connect with foot controller socket in main body.                                |
| Clamping nut                      | Required to fix the plug of the two-button foot controller when attached to the main body. |

Figure 6. Two-button for controller (inside view)

### ③ Micromotor

Figure 7. Micromotor (general view)



Table 5. Description of the micromotor


| Name                           | Description                                                                                                      |
|--------------------------------|------------------------------------------------------------------------------------------------------------------|
| Handpiece Adapter              | A part to connect handpiece                                                                                      |
| Micromotor                     | A BLDC motor to give torque to contra angle handpiece. A rpm and torque are controlled by signal from main body. |
| Micromotor cable and connector | A part to connect with micromotor in main body.                                                                  |

Figure 8. Micromotor (inside view)

④ AC power cable, sterilization plug, handpiece rack, irrigation stanchion, irrigation tube holder

Table 6. AC power cable, sterilization plug, handpiece rack, irrigation stanchion, irrigation tube holder description

| Name                   | Photo                                                                                | Description                                                                                                      |
|------------------------|--------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------|
| AC power cable         |     | Designed to transmit electricity to the main body and is connected to the socket on its rear panel.              |
| Sterilization plug     |    | Fits onto the handpiece adapter during sterilization. Designed to protect the micromotor from steam penetration. |
| Handpiece rack         |   | Stand for positioning the handpiece of the micromotor between use                                                |
| Irrigation Stanchion   |  | Designed to fix the package with the irrigation solution.                                                        |
| Irrigation tube holder |   | Designed to hold the irrigation tubes, installed on the micromotor cable and on the irrigation tube.             |

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- **Dry your hands completely before connecting the power cable.**
- **This product has the free voltage function (AC 100-240 V~, 50-60 Hz).**
- **To avoid risk of electric shock, this equipment must only be connected to supply mains with protective earth.**

#### 8.2.2 Installation of the irrigation system

- Insert the irrigation stanchion in the receptacle on the main body.
- Insert the irrigation tube:
  - Rotate the OPEN/CLOSE lever to "OPEN" position to open the cover for tube insertion.
  - Insert the irrigation tube into the adapter by using both hands. Place the irrigation tube right in the middle of the roller so that it is not kinked or crushed. Turn the OPEN/CLOSE lever to CLOSE after the irrigation tube is placed properly in the tube guide.
- Connect the handpiece with the irrigation tube. More detailed information can be referred to the user manual of handpiece supplier.

#### 8.2.3 Connecting the components

##### Micromotor connection and disconnection

- Locate the marking [↔] of the connector joint upward and push it into the micromotor socket to connect. It will make a "click" sound at the joint when correctly connected. When disconnecting, grab and pull the connector joint. It is imperative to turn off the power supply to the main unit when connecting and disconnecting the plugs.


LCD will show this error message if the micromotor is disconnected or not connected well. Check the connection for normal usage.

### Connecting and disconnecting of two-button foot controller

- To connect, insert the connector of the two-button foot controller into the socket of the two-button foot controller on the back of the main body and turn the clamping nut clockwise.
- To disconnect, turn the clamping nut counterclockwise, grasp the foot controller connector and pull it towards you.

LCD will show this error message if the two-button foot connector is disconnected or not connected well. Check the connection for normal usage.

### 8.2.4 LCD screen application procedure

#### Turn on the unit

Power switch (I/O) can be found in the back side of the product. Press "I".

### Control panel functions

Table 7. Functions and description.

| Button | Description of Function                                                           |
|--------|-----------------------------------------------------------------------------------|
|        | Move the user program ("STEP") up                                                 |
|        | Move the user program ("STEP") down                                               |
|        | Activate the data mode to initiate program modification and save the changed data |
|        | Increase the data cursor is on                                                    |
|        | Decrease the data cursor is on                                                    |

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Table 8. Program (Step) Information\*

| STEP | Ratio | Process  | rpm  | Torque | Pump | Direction |
|------|-------|----------|------|--------|------|-----------|
| 1    | 20:1  | Drilling | 1000 | 35     |      |           |

|   |      |                                      |     |    |     |  |
|---|------|--------------------------------------|-----|----|-----|--|
| 2 | 20:1 | Countersinking                       | 300 | 35 | ●●● |  |
| 3 | 20:1 | Implant Insertion or Tapping         | 20  | 35 | ○●○ |  |
| 4 | 20:1 | Cover screw & Healing cap connection | 100 | 5  | ●●● |  |
| 5 | 20:1 | Removal                              | 20  | 55 | ○●○ |  |
| 6 | 20:1 | Spare                                | 20  | 5  | ○●○ |  |


\* The above table information is not the standard guidance for all implant systems.

### Change STEP

- Program button will move the steps upward (STEP 1-2-3-4-5-6)  
- If you press the button continuously, steps will rotate.  
(The Program move key in the Foot controller unit will make the same change.)
- Program button will move the steps downward (STEP 6-5-4-3-2-1)  
- If you press the button continuously, steps will rotate.

### Change rotation speed

- Press the SELECT button to activate the cursor
- Move the cursor to the rpm change sector by pressing the SELECT button repeatedly.
- Press or button set the target rpm value.
- Press SELECT to save the rpm value displayed on the LCD display. (It will save the value automatically if the value is not changed in 5 seconds.)

The rpm will increase/decrease by the increment value specified in Table 9.

**Table 9. Increment value steps.**


| Setting rpm | Increment rpm |
|-------------|---------------|
| 15-50       | 5             |
| 60-100      | 10            |
| 200-2000    | 100           |

- Rotation speed may become fluctuated ( $\pm 10\%$ ) in the speed below 25 rpm.


### Torque change

- Press the SELECT  button to activate the cursor.
- Move the cursor to the torque change sector by pressing the SELECT button repeatedly.
- Press  or  button to set target torque value.
- Press SELECT to save the torque value displayed on the LCD Display. (It will save the value automatically if the value is not changed in 5 seconds.)
- The torque will increase/decrease by the increment value of 5 N•cm.
- Torque may become fluctuated ( $\pm 10\%$ ) in the torque below 5 N•cm.

### Change direction


- Press the SELECT  button to activate the cursor.
- Move the cursor to the rotating direction change sector by pressing the SELECT button repeatedly.
- Press  or  button to change direction.
- Press SELECT to save the direction displayed on the LCD display. (It will save the value automatically if the value is not changed in 5 seconds.)

- Long beep sound will be given continuously if you put the rotation as reverse direction by the control panel, and short beep sound by the two-button foot controller.
- Changes to the set values can only be made after the micromotor has come to a complete stop.
- "F" stands for Forward, and "R" stands for Reverse.

### Change the gear ratio

- Press the SELECT  button to activate the cursor.
- Move the cursor to the Gear ratio change sector by pressing the "SELECT" button repeatedly.
- Press  or  button to change the gear ratio.
- Press SELECT to save the gear ratio displayed on the LCD display. (It will save the value automatically if the value is not changed in 5 seconds.)

- The available rpm and torque ranges respect to the gear ratio are shown as in Table 15.


**Table 10. Value ranges.**

| Gear ration | rpm     | Torque    |
|-------------|---------|-----------|
| 20:1        | 20-2000 | 5-70 N•cm |
| 27:1        | 15-1400 |           |
| 32:1        | 15-1200 |           |

### Change the irrigation amount

- Press the SELECT button to activate the cursor.
- Move the cursor to the Irrigation amount change sector by pressing the SELECT button repeatedly.
- Press  or  button to change the irrigation amount.
- Press "SELECT" to save the irrigation amount displayed on the LCD display. (It will save the value automatically if the value is not changed in 5 seconds.)

- Water irrigation can be disabled by pressing the  or  button.
- Irrigation amount can be adjusted in 3 steps. Please refer to the below table for the amount of irrigation in each step.

**Table 11. Irrigation amount values**

| Step marking                                                                       | Irrigation amount(Max)<br>mL/min |
|------------------------------------------------------------------------------------|----------------------------------|
|                                                                                    | OFF                              |
|   | 40 mL                            |
|   | 50 mL                            |
|  | 60 mL                            |

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## 9. Technical Specification

**Table 12. Main technical characteristics**

| Parameters                                        | Specification                                                                                                 |
|---------------------------------------------------|---------------------------------------------------------------------------------------------------------------|
| Input power voltage and frequency                 | 100-240 V ~, 50/60 Hz, 130-160 V·A                                                                            |
| Electric shock protection and applied part        | Class 1, type B                                                                                               |
| Handpiece speed and handpiece torque values range | 20-2,000 rpm                                                                                                  |
| Handpiece torque range                            | 5-70 N·cm<br>(based on gear ratio 20:1)                                                                       |
| Waterproof rate of the two-button foot controller | IPX 1                                                                                                         |
| Irrigation speed                                  | Level 1: 40 ml/min<br>Level 2: 50 ml/min<br>Level 3: 60 ml/min                                                |
| LCD screen                                        | Active display area: 320 x 240 mm.<br>Module size: 139 (H) x 100 (W) x 13.1 (D) mm<br>input (typical): 5 VDC. |
| Mode of operation                                 | Short operating mode.<br>Maximum activation time: 3 min<br>Minimum deactivation time: 10 min                  |

|                      |                      |
|----------------------|----------------------|
| Sound pressure level | Not more than 80 dBA |
|----------------------|----------------------|

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**Table 14. Electromagnetic immunity**

| Immunity test                               | IEC 60601 test level                                           | Conformance level                                              | Electromagnetic environment-guidance                                                                                                                        |
|---------------------------------------------|----------------------------------------------------------------|----------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge (ESD) IEC 61000-4-2 | Discharge by 8 kV direct contact<br>15 kV of Air-gap discharge | Discharge by 8 kV direct contact<br>15 kV of Air-gap discharge | Floors should be made of wood or concrete or have ceramic tiles. When the floor is made of synthetic material, the relative humidity must be at least 30 %. |
| Radiated RF Electromagnetic Field immunity  | 3 V/m<br>80 MHz-2.7 GHz<br>80 % AM at 1 kHz                    | 3 V/m<br>80 MHz-2.7 GHz<br>80 % AM at 1 kHz                    | The ICT motor is suitable to use in professional healthcare environment.                                                                                    |

|                                                                                                     |                                                                                                                                                                 |                                                                                                                                                                 |                                                                                                                                                                                                                                                                                                                                |
|-----------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| IEC 61000-4-3                                                                                       |                                                                                                                                                                 |                                                                                                                                                                 |                                                                                                                                                                                                                                                                                                                                |
| Immunity to Proximity Fields from RF wireless Communications Equipment<br>Table 9 in IEC 60601-1-2  | 28 V/m Max.<br>385-5785 MHz in according to table 9                                                                                                             | 28 V/m Max.<br>385-5785 MHz in according to table 9                                                                                                             | RF communication equipment is used on closer than 30 cm to any part of the ICT motor, including cables specified by Dentium.                                                                                                                                                                                                   |
| Electrical fast transient/burst IEC 61000-4-4                                                       | $\pm 2$ kV for power supply lines<br>$\pm 1$ kV for I/O lines<br>100 kHz repetition frequency                                                                   | $\pm 2$ kV for power supply lines<br>$\pm 1$ kV for I/O lines<br>100 kHz repetition frequency                                                                   | The quality of supplied power should be suitable for general business site or hospital environment.                                                                                                                                                                                                                            |
| Surge IEC 61000-4-5                                                                                 | $\pm 1$ kV line to line<br>$\pm 2$ kV line to earth                                                                                                             | $\pm 1$ kV line to line<br>$\pm 2$ kV line to earth                                                                                                             | The quality of supplied power should be suitable for general business site or hospital environment.                                                                                                                                                                                                                            |
| Immunity to conducted disturbances induced by RF fields IEC 61000-4-6                               | 3 V<br>0.15-80 MHz<br>6 V in ISM bands between 0.15 and 80 MHz<br>80 % AM at 1 kHz<br><br>Power supply line & I/O lines                                         | 3 V<br>0.15-80 MHz<br>6 V in ISM bands between 0.15 and 80 MHz<br>80 % AM at 1 kHz<br><br>Power supply line & I/O lines                                         | The strength of RF field in the frequency range higher than 150 kHz~80 MHz, the strength of the RF field is smaller than 3 V.                                                                                                                                                                                                  |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | 0 % UT: 0.5 cycle<br>At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°<br>0 % UT; 1 cycle and 70 % UT; 25/30 Cycles Single phase: at 0°<br>0 % UT; 250/300 cycle | 0 % UT: 0.5 cycle<br>At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°<br>0 % UT; 1 cycle and 70 % UT; 25/30 Cycles Single phase: at 0°<br>0 % UT; 250/300 cycle | The quality of supplied power should be suitable for general business site or hospital environment.<br>For the user to operate the equipment continuously even the electric power supply is interrupted, it is recommended that the uninterruptable power supply device (UPS) or battery is prepared.                          |
| Magnetic field of supply frequency (50/60 Hz) IEC 61000-4-8                                         | 30 A/m                                                                                                                                                          | 30 A/m                                                                                                                                                          | If image distortion occurs, it may be necessary to position the Model ICT motor image intensifier further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low. |

### 9.3 Protection Distance


- Portable and mobile HF communication equipment should be operated at recommended protection distance from this equipment.
- This equipment is manufactured to operate in the electromagnetic environment where the emitted HF disturbance is controlled.

- RF communication equipment is used not closer than 30 cm to any part of the ICT motor including cables specified by Dentium Co.Ltd. (ICT Branch).

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Table 15. Symbols and designations

| Symbol                                                                              | Description                                                                                |
|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|
|    | Type B applied part                                                                        |
|    | Refer to the instructions for use for important information, namely warnings and cautions. |
|    | Caution! Moving machinery                                                                  |
|    | Fuse                                                                                       |
|    | Foot controller (two-button foot controller)                                               |
|   | Do not dispose of with household waste                                                     |
|  | Refer to operation manual/leaflet                                                          |
|  | Caution!                                                                                   |
|  | Note                                                                                       |
|  | Compliance with European standards                                                         |
|  | Medical device manufacturer                                                                |
|  | Date of manufacture of the medical device                                                  |
|  | Authorized Representative in the EU                                                        |
|  | Catalog number                                                                             |
|  | Serial number                                                                              |
| IPX1                                                                                | Water protection degree                                                                    |

|                                                                                     |                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------|
|    | Autoclaving (up to a maximum temperature of 135 °C)   |
|    | Autoclaving (up to a maximum temperature of 134 °C)   |
|    | Top (correct vertical position of the load)           |
|    | Inadmissibility of shock loading                      |
|    | Limit on the number of tiers in a stack               |
|    | Moisture control                                      |
|   | "On" position (state) of the device part power supply |
|  | "On" position (state) of the device part power supply |
|  | Power plug                                            |
|  | Alternating current                                   |
|  | Direct current                                        |
|  | Class II product                                      |

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### Commissioning conditions

- Cables should be checked for peeling of the insulating cover and stored in a roll.
- Stop further use of the device and contact an authorized representative of the manufacturer for service in case of technical problems, such as noise or overheating during operation.
- Dispose of any part of electronic equipment in an environmentally friendly way.
- It is unacceptable to locate the equipment or any part of it in a place exposed to water or moisture.
- Locate the equipment in a location free from exposure to temperature, ventilation, sunlight, dust, salt, etc.
- Avoid placing the equipment on an inclined surface or subjecting it to shock.
- Do not store equipment near chemicals or gases.

- Disconnect power to the equipment and apply to professional assistance if water or liquid enters the equipment.

**Table 16. Checklist of daily check parameters**

| Description                | Check item                                                                                                                                                                                                                                           |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Main body                  | <ol style="list-style-type: none"> <li>1. Is the power switch for unused equipment in the "O" position?</li> <li>2. Is there any damage to the screen panel and buttons?</li> <li>3. Are colors and data clearly displayed on the screen?</li> </ol> |
| Micromotor                 | <ol style="list-style-type: none"> <li>1. Is the micromotor cable damaged?</li> <li>2. Is the micromotor connector damaged?</li> <li>3. Is there any loose screw connection of the micromotor?</li> </ol>                                            |
| Two-button foot controller | <ol style="list-style-type: none"> <li>1. Is the two-button foot controller damaged?</li> <li>2. Does the two-button foot controller move smoothly when pressed and released?</li> </ol>                                                             |

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**Table 17. Possible reasons of malfunctions and solutions.**

| Malfunction                                          | Reason and solution |                                                                                                                                                                                                                                      |
|------------------------------------------------------|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                      | Reason              | Lack of power supply                                                                                                                                                                                                                 |
| LCD does not turn on                                 | Reason              | <ol style="list-style-type: none"> <li>1. Unplug the power connector from the power outlet, wait 5 seconds and reconnect it.</li> <li>2. Turn the power switch of the main unit to the "ON" position. ("I").</li> </ol>              |
|                                                      | Solution            |                                                                                                                                                                                                                                      |
| Buttons do not work                                  | Reason              | <ol style="list-style-type: none"> <li>1. The device is in operating mode or the previous operations are not completed.</li> <li>2. The assigned value is not activated.</li> </ol>                                                  |
|                                                      | Solution            | <ol style="list-style-type: none"> <li>1. Wait until the operation is complete and try again.</li> <li>2. Press the SELECT button to activate the selected value.</li> </ol>                                                         |
| Handpiece does not work                              | Reason              | <ol style="list-style-type: none"> <li>1. Handpiece is defective or needs cleaning.</li> <li>2. Improper connection of a two-button foot controller or micromotor.</li> </ol>                                                        |
|                                                      | Solution            | <ol style="list-style-type: none"> <li>1. Read the operation manual for the handpiece.</li> <li>2. Disconnect and reconnect connectors or check connections.</li> </ol>                                                              |
| Irrigation pump does not work or is not working well | Reason              | <ol style="list-style-type: none"> <li>1. Irrigation tube is damaged or blocked.</li> <li>2. Irrigation tube stopper is blocked.</li> <li>3. Irrigation mode is not activated. Irrigation tube is not connected properly.</li> </ol> |

|  |          |                                                                                                                                                                                                                                                           |
|--|----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Solution | 1. Replace the tube or clean handpiece tip. See handpiece instruction manual for handpiece cleaning.<br>2. Unlock the irrigation tube stopper.<br>3. Activate the irrigation mode. Disconnect and reconnect the two-button foot controller or micromotor. |
|--|----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

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**Table 18. Operation conditions, transportation and storage**

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| Ambient temperature | 10 °C to 35 °C (operation)<br>-10 °C to 60 °C (transportation and storage)                        |
| Relative humidity   | Up to 80 % (operation)<br>Up to 85 % (transportation and storage)                                 |
| Atmosphere pressure | 800-1060 hPa (operation, transportation and storage)<br>700-1060 hPa (transportation and storage) |

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
**Figure 9. Closure of the handpiece connector with a sterilization plug.**

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### **15. List of international standards that the medical device complies with**

EN ISO 13485, EN ISO 14971, EN 60601-1, IEC 80601-2-60, EN 1640, EN 60601-1-2, EN 60601-1-6, EN 62366, EN 62304, ISO 17665-1, ISO 14457, EN ISO 3964, EN ISO 15223-1, EN ISO 9687, EN 1041, EN ISO 14155, MEDDEV 2.12\_1, rev. 8, MEDDEV 2.7.1 rev. 4.

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