HAEMONETICS®

THE Blood Management Company™

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
ENGLISH

PL HIGH EFFICIENCY LEUCOCYTE REMOVAL FILTER FOR PLATELETS (Re-order No: PL1BE)

INDICATIONS

The PL High Efficiency Leucocyte Removal Filter is indicated for the removal of leucocytes, microaggregates and the reduction of anaphylatoxin C3a from a single unit of platelet concentrate prepared from one unit of whole blood or an equivalent volume of single donor platelet collection. A 175 ml PVC administration bag with needleless access port is attached. Sterilized by irradiation. STERILE and non-pyrogenic fluid pathway.

PRECAUTIONS

Follow instructions carefully. For single use; do not resterilize or reuse. Do not use if packaging is damaged or end protectors are loose or displaced. Filtered platelets should be transfused in accordance with standards and/or standard practices.

CAUTION: This device contains di(2-ethylhexyl) phthalate (DEHP). Precautions should be taken to limit long-term exposure to DEHP in certain groups of patients who, based on animal data, may be at risk of adverse effects on reproductive and developmental processes. These are male newborns, infants and young children, periparturient males and pregnant or nursing women. Medical procedures should not be avoided however, as the benefits outweigh any possible health risks associated with DEHP exposure. Please refer to current literature to make an informed decision.

INSTRUCTIONS FOR USE

(Read and follow carefully)

Close clamp (A). Ensure cap of needleless access port (B) is secure. Activate platelet bag and hang on drip stand.

Remove protective cap from filter set spike and insert into platelet bag using a twisting motion.

Open clamp (A) to prime filter.

NOTE: Ensure filter hangs vertically below platelet bag during priming. When filtration is complete, close clamp (A).

Hang administration bag and administer filtered platelets in usual manner.

NOTES:

1. To maintain traceability, do not remove platelet bag from filter set.
2. A needleless access port (B) is provided for syringe access, if required.
3. If an administration set is used for transfusion, any excess air should be removed from the bag by syringe or other appropriate means before transfusion.
4. This product is free of natural rubber latex.

* See product literature for performance summary.

1 C3a is a product of Complement activation.

Haemonetics® and THE Blood Management Company™ are trademarks of Haemonetics Corporation.

Leaflet Part Number: 141PL1BE3, Last Updated: November 2012.

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EC REP

HAEMONETICS UK LTD
5 Ashley Drive, Bothwell,
Scotland, G71 8BS,
UK

CE 0123

i

2

STERILE

STERILE R

PHT DEHP

PHT DEHP


1


2

3

4

Produced by
HAEMONETICS MANUFACTURING INC
1630 Industrial Park Street,
Covina, CA 91722, USA

Assembled in Mexico

Visit us on the Web at
www.haemonetics.com

ENGLISH

RC HIGH EFFICIENCY LEUCOCYTE REMOVAL FILTER WITH ATTACHED SET FOR BLOOD TRANSFUSION

(Re-order No: RC1KLE)

INDICATIONS

The RC High Efficiency* Leucocyte Removal Filter is indicated for the removal of leucocytes, platelets and microaggregates from a single unit of red cells or whole blood.

Sterilised by irradiation.

STERILE and non-pyrogenic fluid pathway.

PRECAUTIONS

Follow instructions carefully.

For single use; do not resterilise or reuse. Do not use if packaging is damaged or end protectors are loose or displaced.

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(Read and follow carefully)

1. PREPARATION OF SET AND BLOOD BAG

Agitate blood bag and hang on drip stand.

Ensure cap (A) on vented spike is closed.

Locate set clamp (B) directly downstream of drip chamber and close completely.

Insert filter set spike into blood bag using a twisting motion.

2. TO PRIME FILTER AND DRIP CHAMBER

Tilt drip chamber to position shown. Squeeze blood bag and open set clamp (B) to prime filter. Continue squeezing blood bag until filter is primed and desired level is obtained in drip chamber, close set clamp (B).

3. TO TRANSFUSE BLOOD

Return drip chamber to hang vertically, as shown. Prime and connect tubing in usual manner.

Regulate flow with set clamp (B). (20 DROPS = APPROXIMATELY 1 mL).

4. TO RECOVER BLOOD

On completion of filtration open cap (A) on vented spike and allow upstream side of filter to drain.

NOTE:

1. The drip chamber should never be squeezed.
2. Saline priming is not required.
3. The filter should not be flushed with saline after filtration.
4. A pressure cuff may be used to prime filter or maintain flow.
5. This product is free of natural rubber latex.

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Leaflet Part Number: 141RC1KLE2, Last Updated: March 2013.

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ENGLISH - LipiGuard® SB REINUSION FILTER FOR SALVAGED BLOOD (Ri-order No: SB1E) is a low priming volume blood filter designed for the reduction of fat particles, emaphilattic C2A, microaggregates and leucocytes from up to 1 unit of washed or unwashed intra-operative or post-operative salvaged blood. The filter is indicated for the reinfusion of blood derived from the surgical site or post-operative wound drainage.

Ri Only.
Sterilized by ethylene oxide.
Spike and non-syringic fluid pathway.

PRECAUTIONS
Follow instructions carefully. For single use, do not sterilize or reuse. Do not use if packaging is damaged or if sealers are loose or displaced.
Pressure Cuff: Use of this filter with a pressure cuff should comply with the recommendations of said salvage equipment manufacturer's instructions for use.
NOTE: This product is not for use with rubber tubes.

INSTRUCTIONS FOR USE
(Read and follow carefully)

1. Locate administration set clamp directly downstream of drip chamber and close completely. Remove short protection cap from outlet port of filter and gentle cap from spike of administration set. Insert administration set spike firmly into outlet port of filter using a holding motion.



2. Agitate blood bag. Remove protective cap from filter spike and insert spike using a twisting motion into one of ports of administration set. Ensure blood bag remains fully seated.

3. Hold filter vertically above the blood bag. (Fig. 1) Squeeze blood bag with one firm medium firm open clamp maintain pressure until desired blood level is obtained in drip chamber and close clamp. Return filter and drip chamber to upright position (Fig. 2).

4. Hang blood bag in open stand. Prime and connect tubing to usual manner.

DO NOT SQUEEZE BUBBLES CHAMBER.

* C2A is a product of ethylene oxide sterilization.
Haemonetics® The Blood Management Company® is LipiGuard® SB (Ri-ordre No: SB1E) and LipiGuard® SB (Ri-ordre No: SB1E) are trademarks registered in the USA.
Lipid Filter Number 02-0422, Order LipiGuard® October 2002

ENGLISH

LipiGuard® SB REINUSION FILTER FOR SALVAGED BLOOD (Ri-order No: SB1E)

INDICATIONS
The LipiGuard SB Reinfusion Filter is a low priming volume blood filter designed for the reduction of fat particles, emaphilattic C2A, microaggregates and leucocytes from up to 1 unit of washed or unwashed intra-operative or post-operative salvaged blood. The filter is indicated for the reinfusion of blood derived from the surgical site or post-operative wound drainage.

Ri Only.
Sterilized by ethylene oxide.
Spike and non-syringic fluid pathway.

PRECAUTIONS
Follow instructions carefully. For single use, do not sterilize or reuse. Do not use if packaging is damaged or if sealers are loose or displaced.
Pressure Cuff: Use of this filter with a pressure cuff should comply with the recommendations of said salvage equipment manufacturer's instructions for use.
NOTE: This product is not for use with rubber tubes.

INSTRUCTIONS FOR USE

- (Read and follow carefully)
1. Locate administration set clamp directly downstream of drip chamber and close completely. Remove short protection cap from outlet port of filter and gentle cap from spike of administration set. Insert administration set spike firmly into outlet port of filter using a holding motion.
 2. Agitate blood bag. Remove protective cap from filter spike and insert spike using a twisting motion into one of ports of administration set. Ensure blood bag remains fully seated.
 3. Hold filter vertically above the blood bag. (Fig. 1) Squeeze blood bag with one firm medium firm open clamp maintain pressure until desired blood level is obtained in drip chamber and close clamp. Return filter and drip chamber to upright position (Fig. 2).
 4. Hang blood bag in open stand. Prime and connect tubing to usual manner.

DO NOT SQUEEZE BUBBLES CHAMBER.

* C2A is a product of ethylene oxide sterilization.
Haemonetics® The Blood Management Company® is LipiGuard® SB (Ri-ordre No: SB1E) and LipiGuard® SB (Ri-ordre No: SB1E) are trademarks registered in the USA.

Lipid Filter Number 02-0422, Order LipiGuard® October 2002

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