76299

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Statement

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Shandong Diano Bio-tech Co.,Ltd declares that we issued Operation manual as attached.

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
Revision date: 2022.06.20

PRODUCT IDENTIFICATION DATA

Type and designation: *Analyzer immunofluorescent*

Model: *FIC-Q100N*

Serial number: _____

Manufacturing date: _____

Name and address

- manufacturer "Shandong Diano Bio-tech Co., Ltd."
Building C6, Chang'an Industrial Park, Cuishan West Street,
Xin'an Street, Anqiu City, Weifang City, Shandong Province, China.

- seller (supplier) "BIO FOCUS" LLC
109387, Russia, Moscow, Tikhaya str., 28, office 216

Phone fax: +7 (916) 900-65-84, +7 (925) 682-12-42, + 86 10 400-053-69-17 / +86 10 010-828-160-22

Website, e-mail:

general issues: www.biofocus.ru

technical issues: _____

service centre: info@biofocus.ru

Analyzer complies with the current regulatory and technical documentation and is recognized as fit for use

Analyzer meets the requirements of regulatory documents:

GOST R 15.013-2016, GOST ISO 13485-2017, GOST R 50444-2020, GOST IEC 61010-2-101-2013, GOST R IEC 62366-2013, GOST R 50267.0-92, GOST 20010-93, GOST R IEC 61326 -1-2014, GOST R IEC 61326-2-6-2014, GOST R IEC 62304-2013, GOST R ISO 14644-4-2002, GOST R ISO 14644-5-2005, GOST R 52539-2006, GOST R 53079.4 -2008, GOST 27025-86, GOST R 52905-2007, GOST 31891-2012, GOST IEC 61010-1-2014, GOST R ISO 18113-1-2015, GOST R ISO 18113-3-2015


ATTENTION OF CONSUMERS!

Thank you for purchasing an Analyzer manufactured by our company!

We ask you to carefully study this Operation Manual, check the correctness and quality of assembly of the Analyzer and its components, ensuring the operation in accordance with the established requirements.

Do not trust the production of these works to random people, avoid independent unskilled actions!

Remember that if you violate the rules of work, you may lose the right to free warranty repair!

ATTENTION!

Only personnel who have studied this Manual and have experience in operating medical equipment are allowed to service the Analyzer.

The manufacturer reserves the right to make design changes that do not lead to deterioration in the technical characteristics, presentation and usability of the Analyzer.

Please send your comments on the improvement of the design and usability to the address of the manufacturer or supplier.

CONTENT

1	General information	5
2	Purpose and scope.....	5
3	Composition and description of the Analyzer	6
4	Basic security measures	26
5	Transportation	30
6	Placement requirements	31
7	Assembly and preparation for operation	31
8	Instructions for use	34
9	Maintenance and repair	53
10	Specific troubleshooting and corrective actions	54
11	Disposal.....	56
12	Warranty obligations.....	57
13	APPENDIX A	59

1 GENERAL INFORMATION

1.1 This Operation Manual (OM) contains a technical description of the Analyzer Immunofluorescence FIC-Q100N (hereinafter referred to as the Analyzer), rules and instructions for its safe operation and other information that users and personnel who carry out transportation, adjustment, maintenance and repairs need to know.

1.2 The purpose of this manual is to provide all pertinent information required for the transport, commissioning, operation and routine maintenance of the Analyzer.

1.3 This manual has been prepared in accordance with GOST R 2.610-2019, GOST R 2.601-2019, GOST IEC 82079-1-2014, GOST R 56018-2014 / ISO / IEC Guide 37: 2012, GOST R 2.105-2019, Order of the Ministry of Health of the Russian Federation dated January 19, 2017 No. 11n "On approval of the requirements for the content of technical and operational documentation of the manufacturer (manufacturer) of a medical device."

1.4 Terms and definitions - according to GOST R IEC / TO 60788-2009, GOST 31508-2012, GOST R 52002-2003, GOST 18996-80, GOST 18322-2016, GOST 27.002-2015, GOST 25866-83, GOST R 55991.7-2014, GOST R 50267.0-92 and GOST R 2.610-2019.

1.5 This Manual should be considered as an integral part of the Analyzer and should always be available and at the disposal of the installer, user (s) and operator (s) for repair and maintenance.

ATTENTION!

In accordance to the Maintenance file "Quantitative detection Kit of SARS-CoV-2 nucleocapsid antigen in human nasal swabs by fluorescence Immunochromatography" "SARS-CoV-2 N Ag".

2 PURPOSE AND SCOPE

2.1. Purpose: The analyzer is designed to process the "Quantitative detection Kit of SARS-CoV-2 nucleocapsid antigen in human nasal swabs by fluorescence Immunochromatography"

“SARS-CoV-2 N Ag” in medical institutions, as well as in medical and research laboratories for in vitro diagnostics in order to detect the SARS-CoV-2 coronavirus and assess the viral load as a test for a new coronavirus infection (COVID-19).

Functional purpose: auxiliary tool in diagnostics. The use of the product has no population and demographic aspects.

Potential consumers, as well as their professional level: In a medical facility of any profile. Clinical laboratory diagnostics doctor, medical laboratory technician.

2.2 The design of the Analyzer complies with the design documentation and meets the specified requirements in accordance with its intended purpose.

By the type of placement, the Analyzer is classified as desktop, closed.

2.3 Operating mode - automatic, repeated-continuous.

2.4 Analyzer according to the classification of GOST 31508-2012, Order of the Ministry of Health of the Russian Federation No. 4n dated June 06, 2012 "On approval of the nomenclature classification of medical devices" (as amended on September 25, 2014) refers to medical devices of potential risk degree 2a, subgroups 5.1, and according to the possible consequences of failure - to products of class B.

Type of medical device in accordance with the nomenclature classification: 325350.

2.5 The analyzer features high test accuracy and speed, excellent stability and low maintenance costs.

3 COMPOSITION AND DESCRIPTION OF THE ANALYZER

3.1 The main functional parts of the Analyzer are shown in Figures 1 and 2.

The analyzer can operate both independently and as part of a complex with an external personal computer.

3.2 The Analyzer includes slots for ID-cards and test-cassettes (ID-card and test-cassette are not included in the set of the medical device "Analyzer immunofluorescent FIC -Q100N", it is included in the set of the medical device" Kit of reagents "SARS-CoV-2 N Ag".

Control is carried out by means of a touch screen; the output of the results is provided both on the screen and by printing on the built-in thermal printer.


Figure 1 - Analyzer

1. Analyzer - inside view after removing the top cover, Figure 2.




Figure 2

2. Description of the main functional elements and components of the Analyzer.


Built-in printer

Front view


Back view


Back side view

Touch screen


Front top view


Bottom view: Analyzer main board, touch screen, touch screen cable to the main board.

ID-card slot


Front view


View from above

Test-cassette reader


Front view

Side view

I / O board


View from above

power supply unit

Name	DJ-U48S-12
Description	Power Supply
Appointment	Analyzer power supply

<i>Specifications</i>	
Color	Black
Material	Fire retardant material
Purity	No greasy stains. No visible contaminants such as clothing fibers, dust, ash, wood tires, chipboard, metal contamination. No marbling / white or black markings /

	inclusions.
Input voltage	$230 \pm 10\%$ (207-253) V AC, 50 Hz.
Output constant voltage, maximum output current	12 V DC, 4 A.
Dimensions (edit)	$110 \pm 5\text{mm} * 60 \pm 5\text{mm} * 34 \pm 5\text{mm}$.
Weight	270 ± 13 gr.
Manufacturer data	Hengyang Dimension Electronic Co., Ltd.


Appearance

network cable

Name	P24 H05VV-F 3x0.75 1.5 m
Description	Network cable
Appointment	Connecting the power supply to the network

Specifications


Color	Black
Material	PVC, copper
Purity	No greasy stains. No visible contaminants such as clothing fibers, dust, ash, wood tires, chipboard, metal contamination. No marbling / white or black markings / inclusions.
Length	1.5± 0.05 m
Weight	190± 9 G.
Manufacturer data	Zhejiang Jintong Nuclear Cable Co., Ltd.

Appearance

USB connection cable

Name	USB - cable UNITEK USB2.0-RS232
Description	Connection cable
Appointment	Data transfer

<i>Specifications</i>	
Color	Brown
Material	PVC, copper, steel
Purity	No greasy stains. No visible contaminants such as clothing fibers, dust, ash, wood tires, chipboard, metal contamination. No marbling / white or black markings / inclusions.
Length	one.5 m
Weight	74 g
Manufacturer data	Beijing Senyuan Chuangda Technology Co., LTD., China

Appearance

printing paper

Name	Printing paper
-------------	----------------

Description	Roll of paper for printing, 57mm wide
Appointment	Print Test Report

<i>Specifications</i>	
Color	White
Material	Thermal paper
Purity	No greasy stains. No visible contaminants such as clothing fibers, dust, ash, wood tires, chipboard, metal contamination. No marbling / white or black markings / inclusions.
Dimensions (edit)	$57 \pm 5 \text{ mm} * 40 \pm 5 \text{ mm}$
Weight	$57 \pm 3 \text{ G.}$
Manufacturer data	Shenzhen Weihuang Printer Co., LTD
Appearance	

3.3 Description of the principle of the analytical method or the principle of operation of the instrument

The principle of testing is photoelectric, using special fluorescent immunoreagents (test cassettes). Studies are usually carried out in the UV region of the spectrum at the wavelength of maximum absorption for a given group of substances, and the radiation is recorded through a filter that does not transmit excitation rays.

The rectangular principle of fluorescence measurement is commonly used:

Light from a radiation source is passed through filters and focused in a flow cell with rectangular light input and output channels. The emission radiation passes through filters and is measured with a photo detector.

NOTE The mobile phase (solvent) should not absorb light either at the absorption wavelength or at the emission wavelength.

Control is carried out by means of a touch screen; the output of the results is provided both on the screen and by printing on the built-in thermal printer.

Description of the target analyte, scientific validity, principle of action

The new coronaviruses are of the β genus. COVID-19 is an acute respiratory infection. People are susceptible to it. Infection with the SARS-CoV-2 antigen is currently the leading cause of COVID-19; asymptomatic infected people can also be a source of infection. According to the current epidemiological study, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include loss of smell, fatigue, fever, and dry cough. In some cases, nasal congestion, runny nose, sore throat, myalgia, and diarrhea are found.

At an early stage of the disease, the total number of leukocytes in the peripheral blood is normal or reduced, and the number of lymphocytes is reduced. In some patients, there was an increase in the level of liver enzymes, muscle enzymes and myoglobin. CRP and erythrocyte sedimentation rate are increased in most patients, while procalcitonin is normal. In severe cases, the amount of D-dimer increases, and peripheral blood lymphocytes gradually decrease. The nucleic acid of the new coronavirus can be found in throat swabs, sputum, lower respiratory secretions, blood and

other samples.

The purpose of the Analyzer Immunofluorescence FIC-Q100N is to determine the amount of the nucleocapsid antigen of the SARS-CoV-2 coronavirus (N-protein) in swabs from the human nose using the classical mechanism : sandwich immunoassay. On the test line (T) there are antibodies "1" to the N-protein of the SARS-CoV-2 virus, antibodies "2" to the N-protein of the SARS-CoV-2 virus, conjugated to fluorescent microspheres fixed on the test line by the sandwich method. Fluorescent microspheres on the test line can be detected by fluorescence immunoassay. The reaction with the sample will also occur on the control line (C), which is coated with goat anti-mouse antibody (GAM).

These products provide quantitative in vitro diagnostics and are applicable to nasal swab specimens.

3.4 The operational and technical characteristics of the Analyzer are shown in Table 3.1.

Table 3.1

Parameter name	Value, norm
1	2
Appearance, installation quality	in accordance with the wiring diagrams, electrical circuit diagrams and reference samples; defects are not allowed
Analyzer supply voltage, V DC:	12
Maximum current, consumed by the Analyzer from the DJ-U48S-12 power supply unit, A:	3
Input voltage of power supply unit DJ-U48S-12, V AC:	$230 \pm 10\%$ (207-253)
Output voltage of power supply unit DJ-U48S-12, V DC:	12
Rated power supply frequency, Hz:	50
Maximum power consumption together with power supply unit, VA:	36
Maximum current consumption of the Analyzer with a power supply unit, A:	0.16
Weight, kg	$1.90 \pm 20\%$
Overall dimensions (length, width, height), mm	$(270 \times 270 \times 120) \pm 5\%$
Touch screen size (diagonal), inches	7
Time of continuous work, h:	8 (followed by a break of 0.5 h)
Maximum time of readiness for work after switching on, min:	1

Number of tests identified	not limited
Number of tested samples	singles, series
Interface for communication with a computer	RS232, USB, network
Supported data transfer	LAN, GPRS, COM
Minimum sample volume for research	according to the characteristics of the test cassettes
units of measurement	
Measuring ranges	
Criteria for the interpretation of the results obtained	
Time from loading a sample into the analyzer until the first lab result is obtained	12 ± 2 sec.
Performance (the number of laboratory analysis results ready for validation, obtained in one hour of analyzer operation)	180 ± 10
Runtime without calibration	Not limited
Repeatability of results (by coefficient of variation CV),%, no more	8
Signal-to-noise ratio,%	0.01 (1: 10000)
Memory capacity by number of measurements	10000
The volume of the characteristic curves (calibration curves) of sets of different series simultaneously stored in the analyzer's memory	5
A substance used as a fluorescent label	europium (Eu)
Requirements for constructive compatibility with the computer	Computers running Windows XP, Windows 7, Windows 10 that support IEEE 802.3 Ethernet data transmission and reception
Radio channel and information transmission modes	GPRS
Protection degree of electrical equipment	IP20

i Notes:

1 Permissible dimensional deviations - according to working drawings, GOST 24643-81,

GOST 25346-2013 and GOST 25347-2013 (not worse than grade 14).

2 Radio frequencies comply with the "Regulations on the procedure for considering materials, conducting an examination and making a decision on the assignment (designation) of radio frequencies or radio frequency channels for radio electronic equipment within the allocated radio frequency bands" dated February 26, 2008 No. 08-23-02-001 (as amended . as of December 20, 2011).

3.5 The analyzer is equipped with the software "DENUO-PEOPLE-v1.3.488", version "DENUO-PEOPLE-v1.3.488" ensuring its operation in the stated modes, compatible with operating systems *Windows XP, Windows 7 and Windows 10.*

The software performs the functions:

- control of the set operating modes;
- display of status and parameters of work;
- execution of the specified algorithms;
- diagnostics of equipment operability with the issuance of messages;
- operational management;
- registration of emergency events;
- formation of a database;
- generating reports;
- changes in control parameters, equipment configuration.

To display information, pop-up windows, trends, graphs of changes in measured technological parameters are used.

3.6 Requirements for the degree of automation of testing and optimization of the laboratory process are met in accordance with GOST R 55991.2-2014.

3.7 Operating Environment of the Analyzer

3.7.1 The analyzer is suitable for operation in a UHL climate, location category 4.2 in accordance with GOST 15150-69 and GOST 15543.1-89.

Permissible ambient temperature during operation: from plus 10 to plus 35 °C, relative humidity: not lower than 10% and not higher than 80% at plus 25 °C, atmospheric pressure: from 84.0 to 106.7 kPa.

3.7.2 Altitude above sea level: no more than 2000 m.

3.7.3 Type of atmosphere –I according to GOST 15150-69.


Characteristics of the environment and aerosols contained in it: non-toxic, non-explosive according to GOST 31610.0-2019, GOST 30852.19-2002, 4th hazard class according to GOST 12.1.007-76, not radioactive according to SanPiN 2.6.1.2523-09.

3.7.4 The analyzer is resistant to impacts caused by falling from a height of 50 mm onto a hard surface in accordance with GOST R 50267.0-92.

3.7.5 The admissible angle of inclination to any side during use is no more than 5°.

3.7.6 The analyzer, when transported in the manufacturer's packaging, retains its properties after exposure to transport shock with a vibration frequency of 2 ... 3 Hz, an acceleration of $30 \text{ m} \cdot \text{s}^{-2}$.

3.7.7 The temperature in the room during the operation of the Analyzer should not change by more than 2°C per hour.

ATTENTION! The use of the Analyzer in other conditions is considered to be improper use.

3.8 Ergonomic and aesthetic characteristics of the Analyzer are generally consistent with GOST 30.001-83, GOST 20.39.108-85, GOST R 56274-2014, GOST R EN 614-1-2003. Switches and switches meet the standards of GOST 23000-78.

3.9 The electrical equipment is installed in accordance with the requirements of the design documentation and the electrical schematic diagram according to RDT 25 106.

3.10 Controls of the Analyzer are supplied with inscriptions (symbols) in accordance with GOST 24505-80, GOST RMEK 60073-2000 and GOST R IEC 878-95, indicating the control area.

Coloring and selection of signal colors are made taking into account the instructions of GOST 12.4.026-2015, GOST R 50267.0-92 and Figure 3.

Symbol	Meaning	Symbol	Meaning
	Carefully! Refer to instructions for use		Refer to instructions for use
	Biological risk		In vitro diagnostic medical device

			Manufacturer

Figure 3 - Symbols available on the Analyzer

3.11 The Analyzer has an indication

- power connections;
- indication of the position of the main switch "ON";
- the occurrence of a malfunction;
- data reception;
- readiness to perform operations;
- the end of operations.

The equipment contains instructions in Russian regarding:

- to the main connections or related aspects;
- to the presence of high voltage;

3.12 The design of the Analyzer is made controllable and maintainable in accordance with GOST 23660-79, GOST R 27.605-2013, GOST 26656-85, accessible to inspection and removal of components to be replaced, and (as necessary) their sanitization.

Repairs are carried out in specialized organizations.

3.13 The analyzer is designed so that its electronic part is maximally adapted to the use of automatic assembly.

The analyzer contains a minimum number of elements for setting the electronic part, but not to the detriment of its technical parameters.

3.14 Analyzer surfaces made of non-corrosion-resistant materials are coated in accordance with GOST 9.032-74, GOST 9.031-74, GOST 9.301-86 (class not lower than III) and GOST 9.104-2018 (condition group 4). Metallic and non-metallic inorganic coatings comply with GOST 9.301-86 and GOST 9.303-84.

3.15. Analyzer reliability indicators

3.15.1 The Analyzer's MTBF is 3000 hours.

Failure of the Analyzer is a violation of its operational state associated with the failure of any component, which entailed deviation of the parameters beyond the declared limits, if at the same time, in order to restore the operational state, it is necessary to replace or repair any component.

3.15.2 The average time to restore the working condition when replacing a faulty part (element) with a serviceable one should be 2 hours.

Duration of daily maintenance - 1 hour, monthly - 5 hours, annual - 10 hours.

3.15.3 Average service life of the Analyzer is 5 years.

The criterion of the limiting state is considered to be technical and economic inexpediency restoration of characteristics.

3.15.4 The coefficient of technical utilization of the Analyzer is 0.90.

The degree of readiness is 0.99 in accordance with GOST 27.002-2015.

3.15.5 The components, devices, apparatus and instruments used in the design have at the time of installation in the Analyzer a residual service life not less than its warranty period.

3.15.6 The Analyzer software has:

- noise immunity, that is, it excludes the possibility of looping and getting into deadlock situations, the ability to function correctly in the event of failures, failures of part of the computing facilities, personnel errors, and also provides an automatic restart when the power supply is restored after a shutdown without issuing false signals and control actions;
- speed, that is, the minimum restart time, reactions to external events, the minimum loss of time for data protection;
- adaptability and development prospects, that is, the simplicity of adapting programs to changes or extensions of the user's tasks without impairing other indicators.

3.15.7 Test cassettes - single use.

3.16 The design of the Analyzer is suitable for disinfection, ensures the drainage of washing liquids; the number of grooves and protrusions that facilitate the accumulation of dirt or body fluids is minimal.

3.17 The analyzer remains operational and electrical safety when biological fluids get on its surface.

The materials are compatible with them.

3.18 The design of the analyzer's components prevents the possibility of their incorrect assembly.


3.19 Patent clearance for manufacturing methods and design solutions is ensured in relation to the Russian Federation and countries where the Analyzer can be supplied, as well as the transfer of technical, information and other documentation.

Joint intellectual property with commercial potential obtained in the course of work must be protected in accordance with Russian and international legislation.

3.20 The Analyzer case has a smooth, uniform, two-tone outer surface, free from burrs, cracks, stains, dirt, scratches, traces of corrosion, pits, dents and chips, in accordance with approved reference standards.

3.21 The materials used are approved for use in medical equipment by the bodies and institutions of Roszdravnadzor.

Test cassettes meet the general standards of GOST R 51088-2013.

ATTENTION!

It is not allowed to remove identification plates from the Analyzer: they bear the serial number and other useful information.

Changes and rebuilding of the Analyzer by the user, not agreed with the manufacturer, are unacceptable.

4 BASIC SAFETY MEASURES

4.1 The analyzer meets the safety requirements when used for the purposes specified in this Manual.

The design ensures the reliability and safety of its operation during the established service life and provides for the possibility of maintenance, sanitary cleaning, repair and operational control.

4.2 During the preparation for operation and operation of the Analyzer, the safety, occupational safety and environmental protection measures specified in this Manual and the Operating Instructions of "Quantitative detection Kit of SARS-CoV-2 nucleocapsid antigen in human nasal swabs by fluorescence Immunochromatography" "SARS-CoV-2 N Ag".

4.3 Safety instructions.

4.3.1 The analyzer meets the standards of constructive and electrical safety in accordance with GOST IEC60065-2013, GOST R ISO / TO 16142-2008, GOST R 50267.0-92 (for class I, type B), GOST R 52770-2016, GOST IEC 61010-1-2014, GOST IEC 61010-2-081-2013, GOST IEC 61010-2-101-2013, GOST IEC 60950-1-2014, GOST 12.2.091- 2012 (Class I, Pollution Degree I, Installation Category II).

4.3.2. Mains cable length : 1 ± 0.05 m.

4.3.3 The temperature of the parts of the Analyzer, accessible for touching without using a tool, during normal operation and at an ambient temperature from plus 10 to plus 35 °C, does not exceed plus 41 °C.

4.3.4 Sanitary and hygienic requirements for the Analyzer are provided in accordance with the "Unified Sanitary and Epidemiological and Hygienic Requirements for Goods Subject to Sanitary and Epidemiological Supervision (Control)" dated May 28, 2010 No. 299, chapter II, section 18. The Analyzer's materials do not provide irritating, toxic and sensitizing effects; no potentially hazardous radiation.

4.3.5 The design of the Analyzer eliminates the possibility of contact with waste and the possibility of operator infection, and also prevents biological contamination.

Note - During use, there is no contact of reagents and biological material with human skin or other organs, so there is no need to carry out toxicological and sanitary tests.

4.3.6 The analyzer meets the requirements of electromagnetic compatibility in accordance with GOST R IEC 61326-1-2014, GOST R IEC 61326-2-6-2014, GOST R IEC 60601-1-2-2014.

The amount of conducted or radiated radio interference generated by the Analyzer does not exceed the values specified GOST 12.1.006-84.

Electromagnetic environment

- Before starting work, it is necessary to assess the risks associated with the influence of electromagnetic fields on the operation of the product.
- Do not use the Analyzer near strong electromagnetic fields (for example, shielded intentional sources of radio frequency radiation), as they may interfere with its normal functioning.
- Application in the vicinity of mobile RF environments. When the test is switched on and when the test is carried out, communication may affect the device, and hence the result itself.

Electromagnetic compatibility

- The analyzer complies with the requirements for immunity and electromagnetic emission given in GOST R IEC 61326-1-2014 and GOST R IEC 61326-2-6-2014.
- The analyzer was designed and tested in accordance with CISPR 11 (class A). In residential areas, it can cause radio interference, in which case you may need to take measures to reduce the interference, such as reorienting, relocating or shielding the location.
- When using external connected devices (test, remote control or monitoring), electromagnetic emission may occur, exceeding the level specified in this standard.

4.3.7 The design of the analyzer's components excludes the possibility of cuts and injuries during its use and repair.

The Analyzer has no accessible metal parts galvanically connected to the power supply.

4.3.8 The design of control electronic devices excludes:

- not envisaged change of operating modes;
- the presence of errors during work;
- the presence of systematic failures.

The electrical equipment of the Analyzer is short-circuit and idle-proof.

4.3.9 Requirements for electrical safety (resistance, insulation strength, grounding and impedance) are provided in accordance with GOST R 50444-2020.

Supply voltage tolerance: $\pm 23V$.

4.3.10 The equivalent sound level during the operation of the Analyzer does not exceed 60 dB in accordance with GOST 12.1.003-2014. Calculation methods - according to GOST 30530-97.

The vibration level during the operation of the Analyzer does not exceed the standards of GOST 12.1.012-2004.

4.3.11 The degree of protection of electrical equipment is IP20 in accordance with GOST 14254-2015.

4.3.12 The Analyzer software has protection against unauthorized configuration by third parties (password protection, encryption, access control) in accordance with GOST R 50739-95 and GOST R 52069.0-2013, anti-virus protection, as well as protection against incorrect use of functions and the input of unacceptable data.

A situation has been eliminated in which it is possible to gain access to the Analyzer's functions bypassing authorization.

4.3.13 Requirements for fire safety of the Analyzer are complied with "Technical Regulations on Fire Safety Requirements" TR RF 005/2008 (Federal Law of the Russian Federation No. 123-FZ of 22.07.2008).

The probability of a fire in (from) the Analyzer is 10^{-6} per unit per year.

4.3.14 The electrical circuit excludes the possibility of spontaneous switching on, off and changing the Analyzer operation modes; complete or partial interruption of power supply and subsequent restoration does not lead to hazardous situations, including spontaneous start-up.

4.4 When using the Analyzer, the following should be provided:

- periodic control over it;
- timely implementation of the regulations for prevention and repair;

- compliance with the requirements for safe placement.

Avoid moisture penetration, as well as shock and strong vibrations by the Analyzer during its transportation.

4.5 In addition to the requirements of this Manual, during the operation of the Analyzer, it is necessary to comply with the requirements of scientific and technical documents of factories - manufacturers of electrical equipment and devices included in it.

4.6 The Analyzer must be repaired by the manufacturer or an authorized organization. Self-elimination of malfunctions (except those regulated in this Manual) and performance of repair and adjustment work is not allowed. Modification of the Analyzer outside the factory is prohibited.

4.7 It is prohibited to use the Analyzer if it does not comply with the requirements of regulatory and technical documentation.

4.8 Do not use third-party reagent kits.

WARNING!

In case of violation of the norms and rules of operation, the requirements of safety measures established in this Manual, even if the violation was single and related only to one of the established requirements, rules and regulations, the manufacturer and seller, regardless of the terms of purchase and the duration of operation of the Analyzer, did not bear any responsibility for its quality and technical condition, as well as for any consequences that occurred during installation and / or during preparation for operation and / or during the operation of the Analyzer, including those that caused damage to the health and life of people, damage environment and human habitat. In case of violation of the manufacturer's rules for using the equipment, the protection applied in this equipment may deteriorate.

DO NOT use the Analyzer:

- *upon detection of mechanical and other damage to components, electrical wiring and switching devices;*
- *in case of damage to structural elements of the equipment.*

IT IS FORBIDDEN to carry out adjustment, repair and similar work when the electrical

equipment is on.

ATTENTION! Failure to comply with these and other safety measures and precautions specified in this Manual may endanger the life and health of the user, and cause accidents.

5 TRANSPORTATION

5.1 Requirements for transportation - according to GOST 23216-78 and GOST R 51908-2002.

Before transporting the Analyzer, make sure that the package is complete and that there are no irregularities in the packaging.

5.2 Component parts of the Analyzer should be protected from movement by gaskets and shock absorbers made of corrugated cardboard.

Replaceable, spare parts and parts are packed in plastic bags.

5.3 The Analyzer must be transported by carriers who specialize in the transportation of the relevant goods and have the proper licenses (permits) and experience in transportation.

5.4 The Analyzer is shipped in a shipping container by any covered vehicle with protection against the effects of atmospheric precipitation and fog, including the following:

- direct transportation by road for a distance of up to 1,000 km on roads with asphalt or concrete surfaces (roads of the first category) without speed limits or at a speed of up to 40 km / h for a distance of up to 250 km on cobblestone or unpaved roads (roads of the second and third categories);

- multimodal transport by rail, air (in heated pressurized compartments), river transport, in combination with each other, and road transport, as well as sea transport.

5.5 The conditions of transportation of the Analyzer in terms of the impact of climatic factors must correspond to group 5 (ОЖ4) in accordance with GOST 15150-69, in terms of the impact of mechanical factors - to group Ж in accordance with GOST 23216-78.

5.6 Loading and unloading of the Analyzer should be carried out in accordance with GOST 12.3.009-76.

The placement and fastening of containers in vehicles should ensure its stable position, exclude the possibility of displacement and impacts against each other, against the walls of the ship-

ping container and vehicles.

5.7 It is forbidden to stack boxes during transportation or place foreign cargo on them.

5.8 The Analyzer should be dispatched to the Far North and equivalent areas in accordance with GOST 15846-2002.

6 PLACEMENT REQUIREMENTS

6.1 Placement of the Analyzer should be carried out taking into account the technical data and safety parameters specified in the technical documentation.

Storage conditions - in group 2 (C) according to GOST 15150-69.

6.2 Storage of the Analyzer is carried out in closed, well-ventilated rooms in conditions excluding moisture condensation, exposure to direct sunlight and corrosive environments (dust, acid and alkali vapours, as well as corrosive gases).

Storage temperature should be within the range from -10 °C to 40 °C, relative humidity: from 20 to 90%, atmospheric pressure: 86... 106 kPa.

Transportation can be performed by any kind of covered transport without moisture, shock and strong vibration at temperatures from -10 °C to 40 °C.

6.3 The distance from the packaging placed on the racks to any object (including the walls and floor of the store) should be at least 0.1 m, the distance to the nearest heater should be at least 1.0 m.

7 ASSEMBLY AND PREPARATION FOR OPERATION

7.1 After unpacking the Analyzer should be placed on the desktop.

Note - It is recommended not to discard *an empty transport box that may be required for future transportation, storage, etc.*

7.2 Next, you should make sure that the Analyzer is complete.

The delivery should include:

- analyzer, model FIC-Q100N (Shandong Diano Bio-tech Co. Ltd., China) - 1 PC;
- power supply unit, model DJ-U48S-12 (Hengyang Dimension Electronic Com. Ltd., China) - 1 pc.;
- network cable, model P24 H05VV-F, Zhejiang Jintong Nuclear Cable Co., Ltd., China - 1

pc .;

- paper for printing, Shenzhen Weihuang Printer Co., LTD., China - 1 pc .;
- cable for USB connection, model UNITEK USB2.0-RS232, Beijing Senyuan Chuangda Technology Co., LTD., China - 1 pc .;
- operation manual - 1 pc.

7.3 The list of reagents and consumables required for measurements is determined by the delivery specification.

7.4 Marking on the package must comply with GOST R ISO 18113-1-2015, GOST R ISO 18113-3-2015, GOST R ISO 15223-1-2020.

7.5 When assembling the Analyzer in order to ensure safety measures, it is prohibited to:

- get down to work without studying the operational documentation;
- start working without checking the technical condition of the Analyzer's elements, its electrical cables and control devices;
- use elements, components, tools and devices not provided for in the design.

7.6 The Analyzer should be installed in winter only after keeping it unpacked in a closed heated room for at least 4 hours.

7.7 The Analyzer is wired according to the electrical schematic diagram.

The safety of the installation of electrical equipment and components must be ensured by meeting the requirements of GOST 12.3.019-80.

7.8 The following requirements are imposed on the placement of the Analyzer:

- it is prohibited to install it at an inclination of more than 5°;
- the mains socket used to connect the mains plug must be in an accessible place to be able to quickly disconnect the Analyzer from the mains;
- the device must be located at a distance of at least 1 m from heating devices;
- it is prohibited to turn on the Analyzer without protective grounding;
- it is not recommended to place or throw heavy objects on the power cable and damage it in any way;
- it is not recommended to place the Analyzer in the following places:
 - with mechanical vibration;
 - with humidity below 20% or above 80% at plus 25 °C;
 - with direct exposure to sunlight;
 - with temperatures below plus 10 or above plus 35 °C;

- do not use the Analyzer in the immediate vicinity of flammable materials;

7.9 The room should be free from aggressive environments, the mass concentration of dust in the air should not exceed $0.75 \text{ mg} / \text{m}^2$, the electric component of the electromagnetic field should not exceed $3 \text{ V} / \text{m}$ in the frequency range from 150 kHz to 300 MHz.

There should be no strong sources of vibration and electromagnetic field around.

7.10 The analyzer should be located on a stable and level surface, free from heavy dust, direct sunlight and corrosive gases.

The surface must be able to support at least 2.5 kg.

7.11 The analyzer should be located in a well-ventilated area with 10 cm of free space around it to ensure trouble-free operation and maintenance.

7.12 When turning on for the first time, connect the power supply to the power connector on the back of the Analyzer, then press the button on the left side of the Analyzer (Figure 3).



Figure 3

WARNING:

The analyzer is grounded via the power supply!

To ensure correct operation of the instrument, the desktop must be connected to the laboratory's grounding system.

7.13 The analyzer must be fully assembled, installed and connected to the mains in accordance with GOST 32144-2013 and GOST R 50571.28-2006 in the following order:

- check the integrity of electrical elements and installation of control devices, the absence of twists in the cables; mount electrical equipment;
- connect all units of electrical equipment to the ground loop;
- check the absence of violations of the insulation of the electrical wiring;
- connect the equipment to the AC mains; when connecting, check the correspondence of the voltage of the network and electrical equipment;
- check the equipment by turning it on.

ATTENTION!

Incorrect or carelessly performed assembly, adjustment, preparation of samples and reagents can lead to complete destruction of the Analyzer or to the fact that it becomes unsafe for the user!

8 INSTRUCTIONS FOR USE

8.1. Before starting operation, control of the correct assembly of individual equipment and installation of test cassettes in accordance with the documentation for them should be ensured.

Make sure by visual inspection that the Analyzer is in good working order.
condition.

8.2 Operation should be carried out in accordance with SP 2.1.3678-20 and GOST ISO / IEC 17025-2019 in clean rooms in accordance with GOST R ISO 14644-4-2002, GOST R ISO 14644-5-2005 and GOST R 52539-2006.

When working with the Analyzer, one should follow the instructions of GOST R 53079.4-2008, GOST 27025-86, GOST R 52905-2007, GOST 31891-2012 and "Rules for the device, safety precautions, industrial sanitation, anti-epidemic regime and personal hygiene when working in laboratories (departments, departments) of sanitary and epidemiological institutions of the system of the Ministry of Health of the USSR ", Moscow, 1981

8.3 During laboratory work, use rubber gloves in accordance with GOST 20010-93, goggles and special laboratory clothing (gown, suit or overalls).

8.4 The analyzer is suitable for operation only with "Quantitative detection Kit of SARS-CoV-2 nucleocapsid antigen in human nasal swabs by fluorescence Immunochromatography"

“SARS-CoV-2 N Ag”, manufactured by Shandong Diano Bio-tech Co. Ltd., China.

8.5 Control of biocontamination in the laboratory should be carried out in accordance with GOST ISO 14698-1-2005 and GOST ISO 14698-2-2005.

ATTENTION!

The analyzer uses chemicals and other hazardous (caustic, irritating and hazardous) chemicals that can be hazardous to health.

The device is designed to handle potentially infectious samples (nasal secretions) that may cause infection / contamination!

8.6 Connecting the analyzer to the electrical network and preparing for operation must be carried out in accordance with the operating manual 26.60.12.004.15580960 OM (section 7).

8.7 Analysis must be carried out on a flat surface with sufficient space around the instrument. It is not allowed to move it during the analysis.

Exposure to dust during analysis is not allowed.

8.8 After turning on the power, the Analyzer goes through initialization and loading. Upon completion of a successful download, it displays its interface (Figure 6).

To enter the service menu, use the password: 123456.


Figure 6

8.9 The Analyzer is operated by pressing the fingers on the touch screen.

8.10 Configuration

8.10.1 Adjustment is carried out in accordance with practical needs before starting work with the Analyzer. After that, you can start testing individual samples, group testing, project management, working with archive records.

8.10.2 To enter the settings mode, press the "System settings" button to display the second-level menu after entering the interface: setting functions, tools, setting the time, information about the device, resetting the settings to factory defaults.

8.10.3 Function setting

You can enter this interface by clicking "Function Settings". The selection of a setting item (Figure 7) is made by means of the icon behind the corresponding item:


Figure 7

8.10.4 Settings

To enter the next interface, you need to click "Settings".

The main tools are: software update, backlight adjustment, quality control, administration, network configuration, port configuration and network connection (Figure 8).


Figure 8

8.10.5 Software update

To go to the next interface, click "Refresh" to enter the next interface (Figure 9).

It is possible to choose to update from the network or USB disk.

When you click Update, you will be given the option to copy the software update package to a USB drive and insert it into the USB port on the back of the analyzer. The analyzer will receive the new version of the program from the USB disk and update the system.

After the update, the Analyzer will reboot.


Figure 9

8.10.6 Adjusting the backlight

To enter the indicated adjustment, press the "Backlight Adjustment" button (Figure 10).


Figure 10

Next - press "Brightness" or "Darker" to adjust the brightness of the backlight.

Then click "Back" to complete the setting and close the dialog box.

8.10.7 Administrator

The Administrator function for commissioning the Analyzer (Figure 11) is available only to professional users.

After entering the password, you can set up auto test, auto test interval, cloud upload and upload sheet, calibration tool, gain adjustment.

8.10.8 Time setting

To enter this interface, you must Press "Time" (Figure 12).


Figure 11


Figure 12

To set the time, click on the field that you want to change, and then the signs "+" or "-" to adjust (Figure 12).

After making changes - click "Save" to confirm the change. Clicking "Cancel" will help cancel the time change.

8.10.9 Analyzer Information

To enter and display the model number and version information of the Analyzer, press "Info" (Figure 13).


Figure 13

8.10.10 Reset to factory settings

When you click "Factory Reset", a dialog box will appear with which you can clear all data and restore factory settings (Figure 14).

After restoring factory settings, the Analyzer automatically restarts.


Figure 14

8.11 Project management

To enter this interface, press "Test tasks". By clicking on the appropriate item, select the information about the study that is suitable for the current check (Figure 15).


Figure 15

- 1 Click Previous Page, Next Page to view the sample test results saved by the analyzer.
- 2 Click "Delete" to delete information about the selected object of verification.

i An example of an ID-card import method:

Step 1: Insert the ID-card into the ID-card slot of the analyzer.

Step 2: In the interface for managing the scanned object, click the "Import ID card" button to complete the import of the scan information.

Additional notes on importing an ID card:

- 1) Each item number and lot number has unique corresponding inspection information data. Re-setting will overwrite the original setting.
- 2) When inserting an ID card into the Analyzer to scan the corresponding barcode and complete one test, it will add the current ID card check item to the list.

8.12 Single sample testing

Step 1: Install the Analyzer correctly and turn it on (Figure 16).


Figure 16

Step 2: Click "Single Sample Test" to switch to the next interface.

Step 3: Select the type of sample to be checked.

Step 4: Confirm the information about the tested sample and the lot number.

The current test item and lot number information are configured in the Project Management section 14.11, which should be consulted when the test item and lot number information need to be changed.

Step 5: After receiving confirmation information - start testing.

Press "Quick Test" on the Analyzer to start testing the test cassette with the sample.

If you press "Standard Test", the countdown will start. The test cassette test will be automatically performed after the countdown ends (Figure 17).


Figure 17

After testing is complete, the result will be displayed (Figure 18).


Figure 18

Step 6: Result of testing.

When the test is completed, the result will be displayed in a list on the screen (Figure 18). By

clicking on the "Curve" it will be possible to display a detailed test result (Figure 19).

Figure 19

8.13 Testing a batch of samples

Batch testing should be performed in the following sequence:

Start incubating -> incubation in progress -> insert the reagent test cassette

-> testing is in progress -> pull out the test cassette

To get started, you need to click "Batch testing" to switch to the desired interface (Figure 20).

Figure 20

Step 1: Add samples.

Click "Add sample" to add data on all samples for testing (Figure 21).


Figure 21

Step 2: start incubation

Add samples to test cassettes. As soon as the application of each of the samples is completed, - press "START INCUBAT." (Figure 22). The current state will change on the screen and the in-

cubation time will begin (Figure 23).


Figure 22


Figure 23

Step 3: start testing

When the incubation time is over, you must insert the test cassette with the sample into the slot of the Analyzer.

Then click "Start test" to start testing (Figure 24).


Figure 24

At the end of the test, the test result is displayed (Figure 25).


Figure 25

When the test is complete, the test result is displayed.

Step 4: Repeat the above steps to test the next sample cassette.

8.14 Measurement records

To initialize this function, press the "History" button on the main menu page (Figure 26).

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Figure 26

To navigate through the list, press the up / down arrows.

Click "Print" to print the selected test result via the built-in printer.

Click "Print" to delete the selected entry.

Click "Delete" to clear the selected entries.

Click "Search" to search for conditions (Figure 27).

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Figure 27

Set time, sample number, analyte. By pressing the "Search" button, find the corresponding entry in the device memory.

8.14 After repair (replacement) or installation, the Analyzer undergoes a partial technical examination in order to establish that the replaced, newly installed or repaired equipment is in good condition, ensuring safe operation.

8.15 Analyzer retrofitted in service, is subject to technical examination in the scope of periodic control with verification of a set of documentation for modernization, visual and measuring control.

8.16 To finish working with the Analyzer, you need to exit the program settings and move the switch to the "Off" position.

8.17 The analyzer, which has worked out its service life, is subjected to examination in order to determine the possibility of extending its safe service life, including:

- visual and measuring control;
- control check of work in all modes;
- determination of the state with the identification of defects, malfunctions, the degree of wear and tear and corrosion;
- testing of safety devices;
- testing of protective grounding (grounding), insulation resistance of electrical networks and electrical equipment, checking the operation of protection in networks with a solidly grounded neutral.

Based on the results of the examination and analysis of the operating conditions, work is being carried out to determine the residual resource of the device and the possibility of extending its safe operation period.

8.18 The owner of the Analyzer must ensure that it is maintained in good condition and safe for operation by arranging proper maintenance. For this purpose, personnel must be appointed who are responsible for organizing and carrying out maintenance and repair work.

The maintenance personnel must know the structure and operation rules of the Analyzer, pass the safety knowledge test and be admitted to independent work in accordance with the established procedure.

ATTENTION!

- disconnect the power cable if the Analyzer is not used for a long period of time;

- when disconnecting the power cable, hold on to the plug and not pull on the cable itself;
- under no circumstances should you disassemble the Analyzer yourself!
- it is necessary to avoid getting water or any other liquid inside the Analyzer case.

8.20 After each working session, it is necessary to check the purity of the analyzer and, if necessary, thoroughly clean it and the corresponding parts with gauze soaked in 0.5% sodium hypochlorous acid solution.

8.21 Taking the Analyzer out of service:

1 Pull out the appliance plug.

2 Thoroughly clean and disinfect the device following the disinfection instructions in Section 9 of this Manual.

3 Pack for planned transportation, for which:

- place the device in the original packing box;
- affix appropriate labels to the outside of the packaging box with sufficient warnings to indicate that

"The product is out of service";

"The product has been in contact with potentially infectious biological substances."

4 Manual transportation of the Analyzer without packing it into its shipping box should be confined to the work area and care should be taken not to damage the instrument.

8.22 Contraindications to the use of the Analyzer are not regulated and are determined by the characteristics of the used reagents (reagents) and test cassettes.

They also impose restrictions on the use.

8.23 The characteristics of the Analyzer are stable under the stated operating conditions under 3.7.

Information on the traceability of the values specified for calibrators or control materials, which is provided by means of available reference methods (methods) of measurement and (or) standards, is determined from the results of clinical trials.

8.24 Conditions necessary for the collection, processing and preparation of samples, data on the stability of the analyzed samples, including the conditions and duration of storage, transportation conditions, as well as imposed restrictions are determined by the characteristics of reagents (reagents) and test cassettes according to 8.23.

8.25 Taking medical decisions based on the results of testing without prior consultation with the attending physician is unacceptable!

In this case, the Code of Medical Ethics of the Russian Federation must be observed.

8.26 Under no circumstances should you disassemble the analyzer yourself!

8.27 The operating time before prophylaxis is from 3 to 6 months, depending on the intensity and conditions of use.

The technical condition is monitored as part of laboratory equipment within the time frame established in a medical and prophylactic institution.

8.28 Disinfectants are used as needed and must be verified to be free of abrasives prior to application.

Instructions for sanitizing (disinfecting) equipment and laboratory premises - SanPiN 1.2.3685-21 "Hygienic standards and requirements for ensuring the safety and (or) harmlessness of environmental factors for humans"

NOTE - Do not use bleach reagents ($\geq 0.5\%$ solution) on any analyzer surface, as oxidizing agents and solvents can damage the analyzer case and touch screen.

Disinfection requirements.

1 Frequency of disinfection.

The FIC-Q100N immunofluorescence analyzer must be disinfected before and after use.

2 Disinfection method.

2.1 Before cleaning and disinfecting the instrument, turn off the power to make sure the power plug is unplugged.

2.2 Wipe the outside of the instrument with a clean, damp cloth. Be careful not to rub the interior or interior surfaces.

2.3 Wipe the surface of the device twice with a clean cloth moistened with 70% (v / v) alcohol solution.

3 List of disinfectants and requirements:

List: 70% ethanol solution.

Requirements: prepare a solution with a concentration of 70% (by volume).

4 Matters requiring attention

4.1 The operator must wear gloves when cleaning and disinfecting.

4.2 In the process of use, it should be borne in mind that alcohol is highly flammable and there should be no open flames.

4.3 People who are allergic to alcohol should be used with caution.

8.29 Specially trained personnel who have been instructed in accordance with GOST 12.0.004-2015, tested for knowledge of the analyzer's design and have electrical safety qualifications of at least III are allowed for maintenance and repairs.

8.30 For providing protection against overcurrent of the mains supply in the event of a failure, the analyzer power supply is equipped with overcurrent protection.

8.31 It is not allowed to make any changes to the device.

The user is responsible for any modifications or changes, as well as for the consequences arising from this.

8.32 During maintenance and repair work, do not smoke and use open fire, leave contaminated rags and other cleaning materials near and on the analyzer equipment.

8.33. There are no conductive circuits and parts in the product that can be exposed to water and other liquids.

8.34. There are no switches on the power supply. The switch is located on the analyzer itself. The switch has an ON / OFF position latching. In the ON position, the switch key is illuminated.

8.35. If the Analyzer is used outside the manufacturer's specifications, the protection level of the equipment may be impaired.

IMPORTANT POINTS when using the Analyzer:

1. Do not place it in places where it is difficult to operate or it is inconvenient to disconnect the device.

2. Do not insert anything into the test cassette slot other than the test cassettes from the SARS-CoV-2 Coronavirus Nucleocapsid Antigen Quantification Kit in Human Nasal Swabs by Fluorescence Immunochromatography SARS-CoV-2 N Ag , manufactured by Shandong Diano Bio-tech Co. Ltd., China

3. If any of the test specimens may become contaminated, protective gloves or other protective measures should be used to avoid skin contact with the opening of the sample addition cassette.

4. Dispose of used test cassettes in accordance with medical waste management regulations to avoid any biohazard.

5. The Analyzer can store the calibration curves related to the product, the batch imported by the manufacturer, and the user test records. When the manufacturer imports via ID-card or transmits via the upper computer any calibration curve directly, the analyzer will automatically save the data.

The analyzer can store calibration data for 5 lots of kits at the same time. Test results during normal user work will be automatically saved in it as a record.

The analyzer can store a maximum of 10,000 records. Saved data can resume automatically after power off without loss. When users select "Factory Reset", the data is deleted.

6. Users using PC software when working with the Analyzer must have the appropriate knowledge: be familiar with such operating systems as Windows XP, Windows7, Windows 10 and with installation, unloading the software installation package into these operating systems.

9 MAINTENANCE AND REPAIR

9.1 Maintenance of the Analyzer as a whole consists in regular inspection of all working units for serviceability and integrity, special attention should be paid to finding faults in electrical connections, damage to protective coatings, housing integrity, electrical wiring insulation, and grounding devices.

Service of the Analyzer is carried out only after it is de-energized.

9.2 The operating time before prophylaxis is from 3 to 6 months, depending on the intensity and conditions of use.

The technical condition is monitored as part of laboratory equipment within the time frame established in a medical and prophylactic institution.

9.3 When using the Analyzer, the external surfaces should be cleaned regularly, but at least once a month.

For cleaning all external parts should use a cloth moistened with 70% ethanol solution. Do not use solvents of any kind (eg.: acetone or gasoline).

List of disinfectants and requirements:

List: 70% ethanol solution.

Requirements: prepare a solution with a concentration of 70% (by volume).

iNOTE -During cleaning, do not wet electrical parts (connectors, etc.). If electrical parts become accidentally wet, dry them thoroughly before turning on the appliance again...

9.4 Replace paper as it runs out.

9.5. Disinfectants are used as needed and must be verified to be free of abrasives prior to application.


9.6. The external surface is cleaned with a cloth moistened with 70% ethanol solution (see above).

Note - Do not use bleach reagents ($\geq 0.5\%$ solution) on any surfaces of the Analyzer, as oxidizing agents and solvents can damage the Analyzer's case and touch screen.

9.7 Parts that are subject to wear or are defective and need to be replaced are requested as indicated:

- customer order numbers;
- type and / or version of the Analyzer;
- its serial number;
- component code numbers;
- component descriptions;
- quantity;
- the name and address of the medical institution, and the delivery address.

ATTENTION!

Replacement of parts with original spare parts guarantees an efficient and long life of the Analyzer.

9.8 Specially trained personnel who have been instructed in accordance with GOST 12.0.004-2015, tested for knowledge of the analyzer's design and have electrical safety qualifications of at least III are allowed for maintenance and repairs.

ATTENTION!

Users can only perform simple routine maintenance operations. In the event of any malfunctions or breakdowns caused by the use of the device, corrective maintenance must be performed by a qualified technician!

The manufacturer accepts no responsibility for unprofessional repair attempts and the resulting consequences!

9.9 During maintenance and repair work, do not smoke and use open fire, leave contaminated rags and other cleaning materials near and on the Analyzer's equipment.

10 SPECIFIC TROUBLESHOOTING AND CORRECTIVE ACTIONS

10.1 Absence of result when pressing the buttons of the Analyzer interface or pressing the buttons of the analyzer indicates a malfunction. It is necessary to restart the instrument by following these steps:

- turn off the power;
- remove the power plug from the socket and plug it in again;
- turn on the power supply.

10.2 Before contacting the service center, it is recommended that you study this section for the subject of determining the possible cause of the problem. A few simple checks and adjustments can eliminate it and lead to normal operation.

10.3 Possible malfunctions and methods of their elimination are shown in Table 10.1.

ATTENTION!

If the above troubleshooting solutions do not help, you need to contact the service department!

Unqualified actions can lead to injury, damage to the Analyzer and void the warranty!

Performing any operation that is expressly prohibited in this Manual, as well as any adjustments, assembly steps not recommended or prohibited in this Manual, will void the warranty!

Table 10.1

Malfunction	Elimination method	Notes (edit)
The analyzer does not turn on	Make sure that the power supply is connected to the analyzer and to the mains, the power button on the analyzer body is pressed. If all these conditions are met, but the analyzer does not turn on, contact the service department.	
Printed information is missing on the printed measurement report	The thermal paper roller is not inserted correctly into the thermal printer. Take out the thermal paper, turn it upside down to the feed roller on the printer cover, put the paper in the slot, close the cover.	
Printer does not print	Make sure the green light on the	

measurement results	top right corner of the thermal printer is on. If it does, make sure there is paper in the printer tray. Install a new roll. If the indicator does not light up, call for service.	
The analyzer issues a measurement error message	The sample cassette is inserted upside down into the analyzer slot. Remove the cassette, insert in the direction of the arrow marked on it, repeat the measurement.	
The touch panel does not respond to finger touches. The Analyzer is unstable.	Turn off the Analyzer, wait 10 seconds, turn it on. If the problem persists, make sure the analyzer is not near high power equipment. If so, move it a few meters away from this equipment, turn on the power. If the problem persists, call for service.	
The LED on the power supply is off.	Replace the power supply. Call for service.	

Notes on common errors:

1. *A user interface error may result in the Analyzer being unable to perform the "test" function. It will not affect the backup and data transfer function. At this stage, you can download the records stored in the analyzer's memory to an external computer, and then transfer the instrument to a representative company for repair.*
2. *The application logic itself is flawed. This error is very unlikely to occur. If this happens, the data from the test performed will be lost. However, the data saved in the Analyzer will not be lost. At this point, turn off the analyzer and restart it for permanent use.*
3. *System or network resources cannot be used, but the test function and data backup function of the analyzer itself are working.*

This may affect the uploading of records to an external computer. When the system or network resources stabilize, this function will automatically resume.

11 DISPOSAL

11.1 In the event of the occurrence of limiting states and a decision that the Analyzer is unsuitable for repair and further operation or the inexpediency of further operation, its equipment must be dismantled and disposed of.

11.2 Prior to disposal, individual devices must be sorted out in order to assess the possibility of further use outside the Analyzer.

11.3 After the end of the service life, if further operation is impossible, the components of the Analyzer after dismantling must be used or disposed of in the prescribed manner in specialized organizations.

The body and other metal parts of the Analyzer are recyclable as secondary resources of ferrous and non-ferrous metals.

11.4 At the end of its service life in medical institutions, the Analyzer is disposed of in accordance with SanPiN 2.1.3684-21 (waste class B), SP 2.1.7.1386-03 and Federal Law "On the Sanitary and Epidemiological Welfare of the Population" No. M 52-FZ dated 30.03. 1999 year

11.5 Waste materials that remain after analysis and cleaning are hazardous, carrying a risk of infection / contamination, therefore, their handling and disposal must comply with the requirements of the current regulatory legal acts.

11.6 The personnel carrying out the disposal must comply with the safety rules for carrying out one or another method of disposal.

12 WARRANTY OBLIGATIONS

12.1 The manufacturer guarantees the compliance of the Analyzer with the requirements of technical and design documentation, subject to the rules of installation, operation, transportation and storage.

12.2 The warranty period of operation is 18 months, calculated from the date of commissioning, but not later than 6 months from the date of shipment to the consumer.

The guaranteed shelf life is 2 years from the date of manufacture.

12.3 During the warranty period, the manufacturer carries out warranty repairs of the Analyzer or out-of-order equipment.

Following repair and / or warranty service, the Analyzer will be returned to the owner in its original condition at the time of shipment.

12.4 The warranty period for components and spare parts is set equal to the warranty period for the Analyzer and expires simultaneously with the expiration of its service life.

12.5 Complaints are submitted during the warranty period, provided that the consumer observes the requirements of the Instructions for use of the Analyzer, with the preparation of a complaint report.

12.6 The warranty does not apply to:

- for scratches, normal wear and tear, wilful damage or damage due to rough handling or non-observance of the operating instructions, improper installation or use of inappropriate voltage;
- damage due to a chemical or electrochemical reaction, rust, corrosion or exposure to water; for damage due to being in an unusual (inappropriate) environment;
- for accidental damage caused by foreign objects or substances, as well as due to cleaning, for repairs carried out by unauthorized service organizations or unauthorized service partners; or repair using non-original spare parts;
- for cables, clamps and electrical connectors;
- at self-repair experiments without the manufacturer's permission;
- in case of damage caused as a result of force majeure (typhoon, fire, flood, hostilities, etc.);
- in cases where the Analyzer is used by unauthorized persons.

12.7 Submission of requests to the Service Center to eliminate the Analyzer's deficiencies is possible only with the simultaneous presentation of a correctly completed warranty card. At the same time, it must clearly indicate information about the Analyzer, the defects in which are to be eliminated (name, serial number), as well as information about the sale of the Analyzer (date of transfer to the buyer, name and address of the seller), certified by the signature and seal (stamp) of the seller, and also buyer's signature.

12.8 The manufacturer reserves the right to make minor changes (modifications) to the Analyzer's design without notice that do not affect the performance of the product as a whole.

APPENDIX A

Analyzer decontamination

A.1 Decontamination (disinfection)

The use of the Analyzer implies the processing of potentially infectious samples (*nasal secretion*) and positive control samples that may cause contamination.

WARNING!

For the above reason, all parts of the device and its accessories must be considered potentially contaminated.

A.2 Decontamination procedure

For decontamination, use a broad spectrum disinfectant (or an equivalent):

70% ethanol solution.

Requirements: prepare a solution with a concentration of 70% (by volume).


ATTENTION!

It is allowed to use disinfectants only on the surfaces of the Analyzer!

When disinfecting, follow these steps:

1. Wear disposable gloves, safety goggles and suitable clothing...
2. Prepare an autoclave bag for all disposable items to be disinfected and label it with autoclave tape.
3. Disconnect the plug to disconnect the appliance from the mains and prevent short circuits....
4. Remove all accessories. All accessories that come with the device must also be disinfected.
5. Gently spray the surface disinfectant. Alternatively, you can wipe the entire instrument with gauze or soft paper soaked in disinfectant solution....
6. After a contact period of at least 10 minutes, repeat the previous steps (1-5) again.
7. After five hours of contact, clean the surfaces of the instrument with a mild detergent or distilled water to remove any disinfectant residue or dirt.
8. Thoroughly dry the Analyzer.

9. Pack the Analyzer and all accessories in the original packaging.
10. Disinfect your hands and wash them with a mild detergent solution.
11. Fill out the disinfection protocol and attach it to the device

A.3 Declaration of decontamination

The declaration must be completed and attached to each device before storage or before being taken to a service center for maintenance or repair.

The declaration must be affixed to the packaging of the device in a conspicuous place.

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(Translation)

(2022) J.J.C.W.J.Z.Zi.No. 234

Applicant: Shandong Diano Bio-tech Co., Ltd., residence: Building C6, Chang'an Industrial Park, Cuishan West Street, Xin'an Street, Anqiu City, Weifang City, Shandong Province, unify the social credit code: 91370784MA3BY89A9U.

Legal Representative: Shengming Yang, male, born on Sept. 2, 1948, the holder of Identity Card No. 110107194809020019.

Item of Notarization: Signature and Seal

This is to certify that Shengming Yang, Legal Representative of Shandong Diano Bio-tech Co., Ltd. came to my office on June 22, 2022 and in front of me affixed his signature and seal of Shandong Diano Bio-tech Co., Ltd. on the Statement in Chinese and English attached hereto.

Notary: Yan Mei

Beijing Jingcheng Notary Public Office

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(Translation)

(2022) J.J.C.W.J.Z.Zi.No. 235

Applicant: Shandong Diano Bio-tech Co., Ltd., residence: Building C6, Chang'an Industrial Park, Cuishan West Street, Xin'an Street, Anqiu City, Weifang City, Shandong Province, unify the social credit code: 91370784MA3BY89A9U.

Legal Representative: Shengming Yang, male, born on Sept. 2, 1948, the holder of Identity Card No. 110107194809020019.

Item of Notarization: Conformity between the translation and the original.

This is to certify that the English version attached hereto is in conformity with the original Chinese version of the Notarial Certificate numbered (2022) J.J.C.W.J.Z.Zi.No. 234

Notary: Yan Mei

Beijing Jingcheng Notary Public Office

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