B. Godwin
26 Nov 2021
BRIAN GODWIN
Notary Public, State of Ohio
My Commission Expires 12-17-2023

«УТВЕРЖДАЮ» / "APPROVE"
ETHICON ENDO-SURGERY, LLC, USA
Associate Director Regulatory Affairs
(должность/position)
Brian K. Rubble
(имя/name)

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INSTRUCTION FOR USE OF OF THE MEDICAL DEVICE

Ultrasonic instruments for generator G11, various modifications,

manufactured by
Ethicon Endo-Surgery, LLC,
475 Calle C, 00969 Guaynabo, Puerto Rico, USA

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2021 year / г.

1. NAME AND DELIVERY SET OF THE MEDICAL DEVICE**1.1 Name of the medical device**

Ultrasonic instruments for generator G11, various modifications (and further in the text: instruments):

1. Harmonic ACE+7 Shears with Advanced Hemostasis for Open and Endoscopic Surgery, shaft length 23 cm composed of:

1.1. Scissors Harmonic ACE + 7 with improved hemostasis function for open and endoscopic operations, barrel length 23 cm - 6 pcs;

1.2. Working part lock gray - 6 pcs.;

1.3. Instructions for use - 1 pc.

2. Harmonic ACE+7 Shears with Advanced Hemostasis for Open and Endoscopic Surgery, shaft length 36 cm, composed of:

2.1. Scissors Harmonic ACE + 7 with improved hemostasis function for open and endoscopic operations, stem length 36 cm - 6 pcs;

2.2. Working part lock gray - 6 pcs.;

2.3. Instructions for use - 1 pc.

3. Harmonic ACE+7 Shears with Advanced Hemostasis for Open and Endoscopic Surgery, shaft length 45 cm, composed of:

3.1. Scissors Harmonic ACE + 7 with the function of improved hemostasis for open and endoscopic operations, stem length 45 cm - 6 pcs;

3.2. Working part lock gray - 6 pcs.;

3.3. Instructions for use - 1 pc.

4. Harmonic ACE+ Shears with Adaptive Tissue Technology for Open and Endoscopic Surgery, shaft length 23 cm, composed of:

4.1. Harmonic ACE + scissors with tissue adaptation technology for open and endoscopic operations, barrel length 23 cm - 6 pcs.;

4.2. Working part lock gray - 6 pcs.;

4.3. Instructions for use - 1 pc.

5. Harmonic ACE+ Shears with Adaptive Tissue Technology for Open and Endoscopic Surgery, shaft length 36 cm, composed of:

5.1. Harmonic ACE + Scissors with Tissue Adaptation Technology for open and endoscopic operations, stem length 36 cm

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1.2 Identification

Harmonic ACE+7 Shears with Advanced Hemostasis for Open and Endoscopic Surgery, shaft length 23 cm	HARH23
Harmonic ACE+7 Shears with Advanced Hemostasis for Open and Endoscopic Surgery, shaft length 36 cm	HARH36
Harmonic ACE+7 Shears with Advanced Hemostasis for Open and Endoscopic Surgery, shaft length 45 cm	HARH45
Harmonic ACE+ Shears with Adaptive Tissue Technology for Open and Endoscopic Surgery, shaft length 23 cm	HAR23
Harmonic ACE+ Shears with Adaptive Tissue Technology for Open and Endoscopic Surgery, shaft length 36 cm	HAR36
Harmonic Focus+ Shears with Adaptive Tissue Technology for Open Surgery	HAR9F
Harmonic Focus+ Long Shears with Adaptive Tissue Technology for Open Surgery	HAR17F
Scissors Harmonic HD 1000i with built-in laparoscopic handle for open and endoscopic operations, stem length 20 cm	SNGHK
Scissors Harmonic HD 1000i with built-in laparoscopic handle for open and endoscopic operations, barrel length 36 cm	SNGCB
Dissection hook blade Harmonic Synergy with a lock of the working part	HARHD20
Curved blade Harmonic Synergy with a lock of the working part	HARHD36
Dissection hook blade Harmonic, diameter 5 mm, length 32 cm	HDH05
Disposable adapter for manual switching, no more than 100 packs	HSA08

2. INFORMATION ABOUT MANUFACTURER MEDICAL DEVICE

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2.1 The manufacturer of the medical device:

«ETHICON ENDO-SURGERY, LLC»

475 Calle C, 00969 Guaynabo, Puerto Rico, USA
tel. +1-513-3371593

2.2 Authorized manufacturer's representative

LLC Johnson & Johnson Krylatskaya 17, Building 2, 121614 Moscow, Russia tel. (495)5807777, fax (495)5807878.

2.3 Manufacturing sites:

1. Ethicon Endo-Surgery, Inc., 4545 Creek Road, Cincinnati, OH 45242, USA;
2. Ethicon Endo-Surgery, Inc., 3801 University Boulevard SE, Albuquerque, NM 87106, USA;
3. Ethicon Endo-Surgery, S.A. de C.V. Planta II, Calle Durango No. 2751, Colonia Lote Bravo, 32575 Ciudad Juarez, Chihuahua, MEXICO;
4. Nipro Healthcare Baja Inc., 3801 University Blvd SE, Albuquerque NM 87106, USA.
5. NPA de Mexico S. de R.L. de C.V., Blvd. Hector Teran Teran #20662-C, Col. Murua Oriente, Tijuana, 22465, Mexico
6. CEA Medical Manufacturing, Inc. d.b.a Nissha Medical Technologies, 1735 Merchants Court Colorado Springs, Colorado 80916 USA

3. DESIGNATION

Designed for dissection of soft tissues in cases where bleeding control and minimization of thermal damage are required.

4. AREA OF APPLICATION

It is used only in medical and preventive health care institutions.

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5. UNDESIRABLE SIDE EFFECTS/RESIDUAL RISKS

Undesirable side effects and risks associated with ultrasonic devices include the potential for bleeding, tissue injury via mechanical or thermal damage, introduction of non-sterile surfaces or pathogen transfer, inflammatory or unintended tissue reaction, electrical shock, foreign body or magnetic resonance incompatibility, and property or environmental damage. Also, unintended harm, extended surgery, or altered surgical approach may result from issues related to device activation, damaged devices, electromagnetic interference, audible noise due to misassembly, misuse of the torque wrench, or an attempt to alter the device.

6. CONTRAINDICATIONS

- The instruments are not indicated for incising bone.
- The instruments are not intended for contraceptive tubal occlusion.
- These instruments are not intended to close the fallopian tubes for contraceptive purposes.

7. USERS

The instruments are intended for use only by highly qualified physicians with sufficient knowledge in the field of surgical techniques.

MEDICAL DEVICE CLASSIFICATION

Depending on the potential risk level of the use in the medical purposes, the medical device "Ultrasonic instruments for generator G11, various modifications" is coded as a Class IIb medical device according to the EC Medical Devices Directive 93/42/EEC.

Type of contact with the human body: devices used for short-term contact with the internal environment and the body tissues. The type of the working part is CF.

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Classification according to suitability for use in an environment with a high oxygen content (according to GOST R IEC 60601-1-2010): the product is not intended for use in an environment with a high oxygen content.

Classification according to the degree of protection against ingress of water according to IEC 60529: IPX0.

Invasiveness: invasive medical device (not implantable).

9. THE POSSIBILITY AND WAYS OF INTEGRATING WITH OTHER MEDICAL DEVICES

The Instruments are designed for use exclusively with:

- the ultrasonic Generator G11 (GEN11) (registration certificate FSZ 2012/11599) software version 2013_1 or later.
- lightweight laparoscopic handle (registration certificate of the FSZ 2011/11397)
- laparoscopic handle for manual activation (registration certificate of the FSZ 2011/11397)
- ultrasonic scalpel "Harmonic" (registration certificate of the FSZ 2012/11532)

10. METHODS AND APPARATUS FOR DISINFECTION AND STERILIZATION CLEANING

Not applicable.

11. INFORMATION IN ACCORDANCE WITH THE DATA OF THE STATE REGISTER OF MEDICINAL PRODUCTS FOR MEDICAL USE ABOUT THOSE CONTAINED IN THE MEDICAL DEVICE

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13. ENVIRONMENT PROTECTION REQUIREMENTS

The instruments do not have a negative impact on people and the environment when used, transported and stored.

• SHELF LIFE

Sterility expiration term is 5 years from date of manufacturing.

15. MANUFACTURER WARRANTY

A storage warranty period - the manufacturer guarantees the preservation of the sterility of the product and its suitability for the intended use, subject to the right transportation and storage within the shelf life as indicated on the product labelling.

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Manufacturer is not liable for any damage or accidents resulting from improper use, human error or external emergencies. Regarding the product quality, please contact an authorized manufacturer's representative.

• MAINTANCE AND REPAIR

The instruments are disposable. Maintenance and repair are not applicable.

17. DISPOSAL

After expiration date of instrument (period of preservation of sterility) and/or in damaged packaging it must be disposed of as medical waste of class A according to SanPiN 2.1.3684-21 "Sanitary-epidemiological requirements for the handling of medical waste" in accordance with the rules of the medical institution.

After use, the product is a potential biological hazard from contamination with blood and other biological materials.

Recycling should also be in accordance with accepted medical practice standards, and relevant local, state, and Federal legislative acts and instructions, as the waste classes B according to SanPiN 2.1.3684-21.

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Ultrasonic instrument for generator G11
*Harmonic Ace+7 Shears with Advanced Hemostasis for Open and
Endoscopic Surgery*
(shaft length 23 cm, 36 cm or 45 cm)

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

Important: This package insert is designed to provide instructions for use of the Harmonic Ace+7 Shears with Advanced Hemostasis for Open and Endoscopic Surgery (shaft length 23 cm, 36 cm or 45 cm). It is not a reference to surgical techniques.

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Harmonic ACE+7 Shears with Advanced Hemostasis for Open and Endoscopic Surgery, shaft length 23 cm	HARH23
Harmonic ACE+7 Shears with Advanced Hemostasis for Open and Endoscopic Surgery, shaft length 36 cm	HARH36
Harmonic ACE+7 Shears with Advanced Hemostasis for Open and Endoscopic Surgery, shaft length 45 cm	HARH45

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Indications

The instruments are indicated for soft tissue incisions when bleeding control and minimal thermal injury are desired. The instruments can be used as an adjunct to or substitute for electrosurgery, lasers and steel scalpels in general, plastic, pediatric, gynecologic, urologic, thoracic, exposure to orthopedic structures (such as spine and joint space), sealing and transection of lymphatic vessels, and other open and endoscopic procedures. The instruments allow for the coagulation of vessels up to and including 7 mm in diameter, using the Advanced Hemostasis hand control button.

Adverse events and complications

Adverse events linked to any surgical procedure or malfunctions are a possibility when and after using this device.

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The following can be considered as an Adverse Event or Malfunction which may occur with use of this device.

Significant Malfunctions:

- Damaged and cracked blade
- Damaged clamp arm
- Melted and detached tissue pad
- Activation issue

Significant Adverse Events

- Burn and injury of Surgeon and patient
- Tissue damage, hemorrhage, leakage\

If an Adverse Event or Malfunction occurs, appropriate techniques or action should be taken.

Description

The Instruments are sterile, single patient use instruments consisting of an ergonomic grip housing assembly with hand control buttons (MIN for minimum power level, MAX for maximum power level, and Advanced Hemostasis for large vessel sealing).

An integrated audible and tactile mechanism in the grip housing indicates full trigger closure. The instruments have a clamp arm and coated curved blade that are designed to work through a 5 mm trocar, through a 5 mm reducer cap in a larger diameter trocar, or through an incision without the use of a trocar. The instrument shafts can be rotated 360° to facilitate visualization and access to targeted tissue. The three dashes on the instrument are intended to represent relative vessel size. The MAX button is typically used for smaller vessels where cutting speed is fastest. The MIN button is typically used in slightly larger vessels and has reduced cutting speed. It is indicated for vessels up to 5 mm in size. The Advanced Hemostasis button is designed for larger vessels and is indicated for vessels up to 7 mm in size. In this mode, cutting speed is further reduced and hemostasis is maximized. The instruments utilize Adaptive Tissue Technology. This provides the generator with the ability to identify and monitor the instrument during

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use, which enables the generator to modulate and adjust its power output as well as provide audible feedback to the user as appropriate.

Each instrument is shipped with one sterile, single-use, disposable torque wrench. Use only the gray torque wrench provided with the instrument. The torque wrench should not be discarded until the completion of the surgical case. Do not attempt to sterilize the disposable torque wrench.

Note: Use of Harmonic torque wrenches other than the one provided may result in damage to the device. The Instruments are designed for use exclusively

with the Generator G11 (GEN11) software version 2013_1 or later. Software revision can be found under

"System Information" in the Generator "Settings" menu. Refer to the Generator Operator's Manual for more information.

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1. Coated Blade/ active jaw (both sides of instrument)
2. Clamp Arm and Tissue Pad
3. Shaft
4. Rotation Knob
5. MAX Hand Control Button (both sides of instrument)
6. MIN Hand Control Button (both sides of instrument)
7. Grip Housing
8. Trigger
9. Hand Piece (not included, registration certificate number FSZ 2011/11397)
10. Torque Wrench
11. Advanced Hemostasis Hand Control Button

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- (3) Shaft – Stainless steel grade UNS S30400
- (1) Blade – Titanium alloy grade Ti 6Al4V ELI
- (1) Blade coating – Fluoropolymer grade Black E2619B
- (2) Tissue Pad – Polytetrafluorethylene grade 7A
- (2) Clamp Arm - Stainless steel grade UNS S30400
- (7) Handle Housing – Polycarbonate grade Calibre 5201-8.

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Storage and Transport Conditions

Temperature: from -22°C to +60°C
Relative Humidity: from 10 to 80% (without condensation)

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Basic specifications

Characteristic	HARH23	HARH36	HARH45
Shaft diameter, mm	5,0	5,0	5,0
Shaft Rotation, °	360°	360°	360°
Shaft length, cm	23	36	45
Targeted vessel max. Diameter, mm	7		
Working temperature range of blades ¹ , °C ²	40 – 300		
Blades' temperature while activation of the improved	40 – 300		

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

hemostasis function control button ¹ , °C ²			
Overall size (A) x (B) x (C), mm	342 x 135 x 32	470x 135 x 32	560 x 135 x 32
Maximum deflection angle of the clamping jaw, °	22°		
Clamp arm grip force (minimum), N	No less than 10		
Length of active jaw (blade), mm	15		
Position of working parts at their closing	"The blade should make contact along the length of the tissue pad when the working parts are closed"		
Sizes for Torque Wrench, mm	50,8 x 68,6 x 25,4		
Mass of Torque Wrench, g	20		
Requirements to the working surfaces and edges of medical instrument	Instrument surfaces shall be free of scratches, nicks, or blemishes		
Requirements to the tightness (impermeability) of working parts for devices	Working parts are sealed		
Requirements to the stability of the working part in biologic fluids	Working parts are resistant to blood		
Minimum width of instrument channel, mm	5,0	5,0	5,0
Mass, g	145	157	160

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Specifications	Value		
	HARH23	HARH36	HARH45
Class of protection against electric shock	Class I		
Type of working part	Type CF		
Operating mode	Continuous duty		
Time of establishment of the operating mode	Max. 10 sec		
Max. output power, W	35		
Max. output voltage, V	150		
Max. loop resistance, Ohm	375		
Maximum current, mA	400		
Blade frequency, Hz	55500		
Max. blade amplitude, μm	85,9	91,2	92,1

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Characteristics of ultrasonic electrosurgical generator GEN11 for connection of instruments

Class of protection against electric shock	Class I
Type of working part	CF type

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Instructions for Use

Verify compatibility of all instruments and accessories prior to using this instrument (refer to Warnings and Precautions).

The Hand Piece is shipped non-sterile. It must be sterilized prior to each use according to the instructions for use supplied with the Hand Piece.

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Assembly

- 1 Using sterile technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field.
- 2 While holding the Hand Piece in a vertical orientation, attach the Hand Piece to the instrument by rotating the instrument onto the Hand Piece in a clockwise rotation as viewed from the distal end of the instrument (finger tight only).

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

- 3 Use the torque wrench (already mounted to the shaft) to tighten the blade onto the Hand Piece. Turn the torque wrench clockwise while holding only the gray Hand Piece until it clicks twice indicating that sufficient torque has been applied to secure the blade.
Note: Do not use any other means than the torque wrench to attach or detach the instrument from the Hand Piece.
Note: Do not torque the instrument by hand without the torque wrench or damage may occur to the Hand Piece.
Note: Hold only the gray Hand Piece and not the instrument handle while applying the torque wrench. (Illustration 2).

- 4 Close the trigger. Remove the torque wrench by sliding it off of the shaft. Do not dispose of the torque wrench until the procedure is completed. The torque wrench is used to remove the instrument from the Hand Piece following the procedure. (Illustration 3) Dispose of the torque wrench only after completing the procedure.
Note: Take care to avoid damage to the blade and clamp arm by closing the trigger while sliding the torque wrench onto or off of the shaft.
Note: Take care to avoid injury from the blade tip while sliding the torque wrench onto or off of the shaft.

Operation

- 1 Connect the assembled Hand Piece and instrument to the ultrasonic generator and turn the ultrasonic generator power on.
Note: MIN is the only mode that can be adjusted.

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- 4 The blade is ultrasonically energized when either the foot switch pedal is depressed or one of the hand control buttons is depressed. Pressing the left foot pedal of the footswitch or the lower hand control button (MIN) on the instrument activates the selected minimum power level. Pressing the right foot pedal of the footswitch or upper hand control button (MAX) on the instrument activates the maximum power level. Pressing the Advanced Hemostasis button

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

on the instrument activates the Advanced Hemostasis. Under default settings, Advanced Hemostasis is hand-activated only.

Note: The "MIN" button can be re-assigned to Advanced Hemostasis by selecting "Convert MIN to Advanced Hemostasis" on the ultrasonic Generator "Settings" screen. Refer to the ultrasonic Generator Operator's

Manual "Settings" section. (Illustration 8) Once "Convert MIN to Advanced Hemostasis" is selected, the Advanced Hemostasis button is disabled.

Note: Energy is not delivered in Advanced Hemostasis mode until the jaws are completely closed.

Caution: If activation is unintentionally stopped while sealing, maintain jaw closure and reactivate.

Note: The ultrasonic generator provides an audible tone to indicate when the instrument blade is active. The ultrasonic generator changes to a second activation tone as Adaptive Tissue Technology regulates the delivery of energy.

Thermal influences such as fluids or minimal to no tissue in the jaws may affect the presence or timing of the tone change.

The tone change does not provide confirmation of tissue effect. When the second tone is heard, the situation should be assessed and the intended surgical action completed, such as gradual application of tension to facilitate transection.

The secondary activation tone change is not a substitute for surgical experience.

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- 5 Close the clamp arm by closing the trigger, and insert the shaft through a trocar or incision.
- 6 The instrument can be used for dissection, grasping, coagulation, and cutting between the blade and clamp arm.
Note: To achieve complete sealing, the trigger should be fully closed and the vessel fully contained between clamp arm and blade of device. An audible and tactile "click" indicates full trigger closure. To achieve full closure of the jaws of the device, squeeze the plastic trigger until you feel it stop against the plastic handle (plastic to plastic) (Illustration 9). If full trigger closure is released prior to or during activation on tissue, an audible and tactile "click" is evident. Increase grip force until full trigger closure is achieved.
Note: Keep clamp arm open when using the inside bottom of the blade for backcutting. (Illustration 10)
WARNING: Do not use Advanced Hemostasis mode for procedures where energy application is desired prior to full closure of the jaws (e.g. solid organs). Energy is not delivered in Advanced Hemostasis mode until the jaws are completely closed.
WARNING: During benchtop testing of vessels >5 mm, the strongest vessel seals were achieved by allowing the Advanced Hemostasis mode to completely transect the targeted vessel.
WARNING: Prolonged usage of Advanced Hemostasis Mode may cause tissue pad damage.

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- 7 Close the clamp arm by closing the trigger and remove the shaft from the trocar or incision.

Shaft Rotation

The instrument's shaft can be rotated 360° to facilitate visualization and access to target tissue when dissecting, grasping, coagulating, and cutting.

Disassembly

- 1 Turn the ultrasonic generator OFF at the power switch.
- 2 Close the clamp arm and slide the torque wrench over the distal end and up the shaft until the wrench aligns with the flats on the shaft. Hold by the Hand Piece only, not the instrument handle, and loosen the instrument by turning the wrench counterclockwise. Continue to loosen by turning the rotation knob manually to completely unscrew the instrument. Do not untorque the instrument by hand without the torque wrench or damage may occur to the Hand Piece.
Note: Do not use any other means than the torque wrench to detach the instrument from the Hand Piece.
Note: Take care to avoid injury from the blade tip while sliding the torque wrench onto or off of the shaft.
- 3 Remove the torque wrench by pulling it straight back over the blade.

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4 Dispose of the instrument in an appropriate container.

Warnings and Precautions

- **ATTENTION!** Product modifications are not allowed! (Applicable for all versions of ultrasonic instruments for the electrosurgical ultrasound generator GEN11)
- Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary in diameter from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- A thorough understanding of the principles and techniques involved in laser, electrosurgical, and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that electrical insulation or grounding is not compromised. Do not immerse instruments in liquid unless the instruments are designed and labeled to be immersed.
- Verify compatibility with generators. Use device only with Ethicon Endo-Surgery Generator G11 (GEN11) software version 2013_1 or later. Software revision can be found under "System Information" in the Generator G11 (GEN11) "Settings" menu. Refer to the Generator G11 (GEN11) Operator's Manual for more information.
- In case of system failure, ensure the availability of the appropriate back up equipment relevant to the specific procedure.
- Audible high-pitched ringing, resonating from the blade or Hand Piece, are an abnormal condition and an indicator that the blade or Hand Piece is not operating properly. The ringing may be an indicator that the Hand

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Piece is beyond its useful life or that the blade has not been attached properly, which may result in abnormally high shaft temperatures and user or patient injury.

- The instruments allow for the coagulation of vessels up to and including 7 mm in diameter, using the Advanced Hemostasis hand control button. Do not attempt to seal vessels in excess of 7 mm in diameter.
- Blood and tissue buildup between the blade and shaft may result in abnormally high temperatures at the distal end of the shaft. To prevent burn injury, remove any visible tissue buildup at the distal end of the shaft.
- As with all energy sources (Electrosurgery, Laser, or Ultrasound), there are concerns about the carcinogenic and infectious potential of the by-products, such as tissue smoke plume and aerosols. Appropriate measures such as protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.
- Do not attempt to bend, sharpen, or otherwise alter the shape of the blade. Doing so may cause blade failure and user or patient injury.
- To avoid user or patient injury in the event that accidental activation occurs, the instrument blade, clamp arm, and distal end of the shaft should not be in contact with the patient, drapes, or flammable materials while not in use.
- During and following activation in tissue, the instrument blade, clamp arm, and distal 7 cm of the shaft may be hot. Avoid unintended contact with tissue, drapes, surgical gowns, at all times.
- Avoid contact with any and all metal or plastic instruments or objects when the instrument is activated. Contact with staples, clips or other instruments while the instrument is activated may result in cracked or broken blades.
- Do not introduce or withdraw the instrument with the jaws open through a trocar sleeve as this may damage the instrument.
- Care should be taken not to apply pressure between the instrument blade and tissue pad without having tissue between them. Clamping the tissue pad against the active blade without tissue on the full length of the blade will result in higher blade, clamp arm and distal shaft temperatures and can result in possible damage to the instrument. If this occurs, there

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

may be an instrument failure, and the generator touchscreen displays a troubleshooting message.

- To avoid user or patient injury, do not activate an electrosurgical device in close proximity to the HARMONIC instruments. The aerosols created by the activation of the HARMONIC instruments in fatty tissue are potentially flammable.

- Keep the clamp arm open when backcutting or while the blade is active without tissue between the blade and tissue pad to avoid damage to the tissue pad and increased blade, clamp arm and distal shaft temperatures.

- The entire exposed blade tip and any exposed blade shaft are active and will cut/coagulate tissue when the instrument blade is activated. Be careful to avoid inadvertent contact between all exposed blade surfaces and surrounding tissue when using the instrument.

- Use only the appropriate Foot Switch, Hand Piece, instruments, and power cord to ensure that they are compatible with the generator.

- After removing the instrument, examine the tissue for hemostasis. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.

- Successful hemostasis may require adjunct measures when HARMONIC instruments are used on solid organs. Due to the difficulty of visualizing internal structures, proceed slowly and do not attempt to transect large masses of tissue in one activation. Avoid the division of large vascular/biliary bundles when using the instrument under these conditions.

- If activation is unintentionally stopped while sealing, maintain jaw closure and reactivate.

- Do not use Advanced Hemostasis mode for procedures where energy application is desired prior to full closure of the jaws (e.g. solid organs). Energy is not delivered in Advanced Hemostasis mode until the jaws are completely closed.

- During benchtop testing of vessels >5 mm, the strongest vessel seals were achieved by allowing the Advanced Hemostasis mode to completely transect the targeted vessel.

- Prolonged usage of Advanced Hemostasis Mode may cause tissue pad damage.

- Products manufactured or distributed by companies not authorized by Ethicon Endo-Surgery may not be compatible with the HARMONIC

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

system. Use of such products may lead to unanticipated results and possible injury to the user or patient.

- If during use the generator displays "Advanced Features Are Not Available In This Device", the following functions of Adaptive Tissue Technology are no longer available: Regulated Energy Delivery, Enhanced Audible Feedback and Advanced Hemostasis. As a consequence, the device vessel sealing indication does not exceed 5 mm.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Incidental and prolonged activation against solid surfaces, such as bone, may result in blade heating and subsequent blade failure, and should be avoided.
- Dispose of all opened instruments whether used or unused. This device is packaged and sterilized for single use only.
- Reuse and improper reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness, or death.
- Reuse and improper reprocessing or resterilization of single-use devices may create a risk of contamination and/or cause infection or cross-infection, including, but not limited to, the transmission of infectious diseases. Contamination may lead to injury, illness, or death.

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

Ethicon
Endo-Surgery | Harmonic

Ultrasonic instrument for generator G11
*Harmonic ACE+ Shears with Adaptive Tissue Technology for Open and
Endoscopic surgery*
(shaft length 23 cm or 36 cm)

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Important: This package insert is designed to provide instructions for use of the Harmonic ACE+ Shears with Adaptive Tissue Technology for Open and Endoscopic surgery (shaft length 23 cm or 36 cm). It is not a reference to surgical techniques.

Instructions for use

This IFU applies to ultrasonic Instruments for ultrasonic electrosurgical generator GEN11, with the following modifications:

Harmonic ACE+ Shears with Advanced Hemostasis for Open and Endoscopic Surgery, shaft length 23 cm	HAR23
Harmonic ACE+ Shears with Advanced Hemostasis for Open and Endoscopic Surgery, shaft length 36 cm	HAR36

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Indications

The Instruments are indicated for soft tissue incisions when bleeding control and minimal thermal injury are desired. The instruments can be used as an adjunct to or substitute for electrosurgery, lasers and steel scalpels in general, plastic, pediatric, gynecologic, urologic, thoracic, exposure to orthopedic structures (such as spine and joint space) and other open and endoscopic procedures.

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•

Device Description

The Instruments are sterile, single patient use instruments consisting of an ergonomic grip housing assembly with hand control buttons (MIN for minimum power level and MAX for maximum power level).

An integrated audible and tactile mechanism in the grip housing indicates full trigger closure. The instruments have a clamp arm and coated curved blade that are designed to work through a 5 mm trocar, through a 5 mm reducer cap in a larger diameter trocar, or through an incision without the use of a trocar. The instrument shafts can be rotated 360° to facilitate

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visualization and access to targeted tissue. The instruments allow for the coagulation of vessels up to and including 5 mm in diameter. Adaptive Tissue Technology provides the ultrasonic generator with the ability to identify and monitor the instrument during use, which enables the ultrasonic generator to modulate and decrease its power output as well as provide audible feedback to the user as appropriate.

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Note: Use of HARMONIC torque wrenches other than the one provided may result in damage to the device.

The instruments are designed for use exclusively with the ultrasonic Generator G11 (GEN11) software version X or later. Software revision can be found under "System Information" in the ultrasonic Generator "Settings" Menu. Refer to the ultrasonic Generator Operator's Manual for more information.

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1. Coated Blade (active jaw);
2. Clamp Arm and Tissue Pad;
3. Shaft;
4. Rotation Knob;
5. MAX Hand Control Button;
6. MIN Hand Control Button;
7. Grip Housing;
8. Trigger;
9. Hand Piece (not included, registration certificate number FSZ 2011/11397);

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- (3) Shaft – Stainless steel grade UNS S30400
- (1) Blade – Titanium alloy grade Ti 6Al4V ELI
- (1) Blade coating – Fluoropolymer grade Black E2619B
- (2) Tissue Pad – Polytetrafluorethylene grade 7A
- (2) Clamp Arm - Stainless steel grade UNS S30400
- (7) Handle Housing – Polycarbonate grade Calibre 5201-8

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Basic specifications

Base specifications	HAR23	HAR36
Shaft diameter, mm	5,0	5,0
Shaft Rotation, °	360°	360°
Shaft length, cm	23	36
Targeted vessel max. Diameter, mm	5	
Working temperature range of blades ¹ , °C ²	40 – 300	
Overall size (A) x (B) x (C), mm	342x 134 x 32	470 x 135 x 32
Maximum deflection angle of the clamping jaw, °	22°	
Clamp arm grip force (minimum), N	no less than 10	

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Lenght of active jaw (blade), mm	15	
Position of working parts at their closing	"The blade should make contact along the length of the tissue pad when the working parts are closed"	
Sizes for Torque Wrench, mm	50,8 x 68,6 x 25,4	
Mass of Torque Wrench, g	20	
Requirements to the working surfaces and edges of medical instrument	Instrument surfaces shall be free of scratches, nicks, or blemishes	
Requirements to the tightness (impermeability) of working parts for devices	Working parts are sealed	
Requirements to the stability of the working part in biologic fluids	Working parts are resistant to blood	
Minimum width of instrument channel, mm	5,0	5,0
Mass, g	145	157

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Specifications	Value	
	HAR23	HAR36
Class of protection against electric shock	Class I	
Type of working part	Type CF	
Operating mode	Continuous duty	
Time of establishment of the operating mode	Max. 10 sec	
Max. output power, W	35	
Max. output voltage, V	150	
Max. loop resistance, Ohm	375	
Maximum current, mA	400	
Blade frequency, Hz	55500	

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Preventive inspection

Before use, be sure to check the tool or any other used medical devices for the absence of any defects or damage, e.g., rough surfaces, nicks or blemishes.

Assembly

- 1 Using sterile technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field.
- 2 While holding the Hand Piece in a vertical orientation, attach the Hand Piece to the instrument by rotating the instrument onto the Hand Piece in a clockwise rotation as viewed from the distal end of the instrument (finger tight only).
- 3 Use the torque wrench (already mounted to the shaft) to tighten the blade onto the Hand Piece. Turn the torque wrench clockwise by holding only the gray Hand Piece until it clicks twice indicating that sufficient torque has been applied to secure the blade.
Note: Do not use any other means than the torque wrench to attach or detach the instrument from the Hand Piece.
Note: Do not torque the instrument by hand without the torque wrench or damage may occur to the Hand Piece.
Note: Hold only the gray Hand Piece and not the instrument handle while applying the torque wrench. (Illustration 2).
- 4 Close the trigger. Remove the torque wrench by sliding it off of the shaft. Do not dispose of the torque wrench until the procedure is completed. The torque wrench is used to remove the instrument from

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

the Hand Piece following the procedure. (Illustration 3) Dispose of the torque wrench only after completing the procedure.

Note: Take care to avoid damage to the blade and clamp arm by closing the trigger while sliding the torque wrench onto or off of the shaft.

Note: Take care to avoid injury from the blade tip while sliding the torque wrench onto or off of the shaft.

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Operation

- 1 Connect the assembled Hand Piece and instrument to the ultrasonic generator and turn the ultrasonic generator power on.

Note: MAX power is set at power level 5 and cannot be adjusted.

- 2 Select the desired minimum power level using the INCREASE and DECREASE buttons on the ultrasonic generator touchscreen.

Note: The recommended minimum starting power level is Level 3 (Illustration 4). For greater tissue cutting speed use a higher ultrasonic generator power level, and for greater coagulation use a lower ultrasonic generator power level. The amount of energy delivered to the tissue and resultant tissue effects are a function of many factors, including the power level selected, blade characteristics, grip force, tissue tension, tissue type, pathology, and surgical technique.

- 3 For optimal performance and to avoid tissue sticking, clean the instrument blade, clamp arm, and distal end of the shaft throughout the procedure by activating the instrument tip in saline. (Illustration 5)

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Note: Do not touch the instrument to metal while activated. (Illustration 6)

Note: Do not clean the blade tip with abrasives. It can be wiped with a moist gauze sponge to remove tissue, if necessary.

If tissue is still visible in the clamp arm, use hemostats to remove residue, taking care not to actuate the Hand Piece. If desired, the instrument may be unplugged. (Illustration 7)

- 4 The blade is ultrasonically energized when either the foot switch pedal is depressed or one of the hand control buttons is depressed. Pressing the left foot pedal of the footswitch or the lower hand control button (MIN) on the instrument activates the selected minimum power level. Pressing the right foot pedal of the footswitch or upper hand control button (MAX) on the instrument activates the maximum power level.

Note: The ultrasonic generator provides an audible tone to indicate when the instrument blade is active. The ultrasonic generator changes to a second activation tone when there is little or no tissue between the instrument tissue pad and blade.

The tone change does not necessarily provide confirmation of tissue effect. Tissue effects are a function of many factors, including the ultrasonic generator power level selected, grip force, tissue tension, tissue type, pathology, and surgical technique, among others. When the second tone is heard, the operative situation should be assessed and the intended action upon the tissue completed accordingly, such as through the gradual application of tension in order to facilitate transection.

Adaptive Tissue Technology determines minimal tissue conditions by monitoring changes in the thermal condition of the blade. As such, thermal influences such as cool fluids or tissue touching the underside or non-clamping surfaces of the blade may affect the presence or timing of the tone change.

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

The activation tone change is not a substitute for surgical experience and should always be interpreted within the context of traditional surgical cues, such as visual or tactile observations, and clinical experience.

The second activation tone can be turned off under "Tone Change ON/OFF" in the ultrasonic Generator "Settings" Menu.

Note: This will deactivate the second activation tone only; this will not affect the Adaptive Tissue Technology's modulation and decrease of power output.

Note: Scratches on the blade may lead to premature blade failure.

Avoid accidental contact with other instruments during use.

Do not use any means other than the torque wrench to attach or detach the instrument from the Hand Piece.

5 Close the clamp arm by closing the trigger, and insert the shaft through a trocar or incision.

6 The instrument can be used for dissection, grasping, coagulation, and cutting between the blade and clamp arm.

Note: Optimum instrument performance is achieved with full trigger closure. An audible and tactile "click" indicates full trigger closure. To achieve full closure of the jaws of the device, squeeze the plastic trigger until you feel it stop against the plastic handle (plastic to plastic) (Illustration 8). If full trigger closure is released prior to or during activation on tissue, an audible and tactile "click" is evident. Increase grip force until full trigger closure is achieved.

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Note: Keep clamp arm open when using the inside bottom of the blade for backcutting. (Illustration 9)

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Shaft Rotation

The instrument's shaft can be rotated 360° to facilitate visualization and access to target tissue when dissecting, grasping, coagulating, and cutting

Disassembly

- 1 Turn the ultrasonic generator OFF at the power switch.
- 2 Close the clamp arm and slide the torque wrench over the distal end and up the shaft until the wrench aligns with the flats on the shaft. Hold by the Hand Piece only, not the instrument handle, and loosen the instrument by turning the wrench counterclockwise. Continue to loosen by turning the rotation knob manually to completely unscrew the instrument. Do not untorque the instrument by hand without the torque wrench or damage may occur to the Hand Piece.
Note: Do not use any other means than the torque wrench to detach the instrument from the Hand Piece.
Note: Take care to avoid injury from the blade tip while sliding the torque wrench onto or off of the shaft.

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3 Remove the torque wrench by pulling it straight back over the blade.

4 Dispose of the instrument in an appropriate container

Warnings and Precautions

- **ATTENTION!** Product modifications are not allowed! (Applicable for all versions of ultrasonic instruments for the electrosurgical ultrasound generator GEN11)
- **Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician.
- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary in diameter from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- A thorough understanding of the principles and techniques involved in laser, electrosurgical, and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that electrical insulation or grounding is not compromised. Do not immerse instruments in liquid unless the instruments are designed and labeled to be immersed.
- Verify compatibility with generators. Use device only with Ethicon Endo-Surgery Generator G11 (GEN11) software version X or later. Software revision can be found under "System Information" in the Generator G11 (GEN11) "Settings" Menu. Refer to the Generator G11 (GEN11) Operator's Manual for more information.
- In case of system failure, ensure the availability of the appropriate back up equipment relevant to the specific procedure.
- Audible high-pitched ringing, resonating from the blade or Hand Piece, are an abnormal condition and an indicator that the blade or Hand Piece is not operating properly. The ringing may be an indicator that the Hand Piece is beyond its useful life or that the blade has not been attached properly,

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which may result in abnormally high shaft temperatures and user or patient injury.

- The instruments allow for the coagulation of vessels up to and including 5 mm in diameter. Do not attempt to seal vessels in excess of 5 mm in diameter.

- Blood and tissue buildup between the blade and shaft may result in abnormally high temperatures at the distal end of the shaft. To prevent burn injury, remove any visible tissue buildup at the distal end of the shaft.

- As with all energy sources (Electrosurgery, Laser, or Ultrasound), there are concerns about the carcinogenic and infectious potential of the by-products, such as tissue smoke plume and aerosols.

2 Appropriate measures such as protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.

- Do not attempt to bend, sharpen, or otherwise alter the shape of the blade. Doing so may cause blade failure and user or patient injury.

- To avoid user or patient injury in the event that accidental activation occurs, the instrument blade, clamp arm, and distal end of the shaft should not be in contact with the patient, drapes, or flammable materials while not in use.

- During and following activation in tissue, the instrument blade, clamp arm, and distal 7 cm of the shaft may be hot. Avoid unintended contact with tissue, drapes, surgical gowns, at all times.

- Avoid contact with any and all metal or plastic instruments or objects when the instrument is activated. Contact with staples, clips or other instruments while the instrument is activated may result in cracked or broken blades.

- Do not introduce or withdraw the instrument with the jaws open through a trocar sleeve as this may damage the instrument.

- Care should be taken not to apply pressure between the instrument blade and tissue pad without having tissue between them. Clamping the tissue pad against the active blade without tissue on the full length of the blade will result in higher blade, clamp arm and distal shaft temperatures and can result in possible damage to the instrument. If this occurs, there may be an instrument failure, and the generator touchscreen displays a troubleshooting message.

- To avoid user or patient injury, do not activate an electrosurgical device in close proximity to the HARMONIC instruments. The aerosols created by

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

the activation of the HARMONIC instruments in fatty tissue are potentially flammable.

- Keep the clamp arm open when backcutting or while the blade is active without tissue between the blade and tissue pad to avoid damage to the tissue pad and increased blade, clamp arm and distal shaft temperatures.
- The entire exposed blade tip and any exposed blade shaft are active and will cut/coagulate tissue when the instrument blade is activated. Be careful to avoid inadvertent contact between all exposed blade surfaces and surrounding tissue when using the instrument.
- Use only the appropriate Foot Switch, Hand Piece, instruments, and power cord to ensure that they are compatible with the generator.
- After removing the instrument, examine the tissue for hemostasis. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.
- Successful hemostasis may require adjunct measures when HARMONIC instruments are used on solid organs. Due to the difficulty of visualizing internal structures, proceed slowly and do not attempt to transect large masses of tissue in one activation. Avoid the division of large vascular/biliary bundles when using the instrument under these conditions.
- Products manufactured or distributed by companies not authorized by Ethicon Endo-Surgery may not be compatible with the HARMONIC system. Use of such products may lead to unanticipated results and possible injury to the user or patient.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Incidental and prolonged activation against solid surfaces, such as bone, may result in blade heating and subsequent blade failure, and should be avoided.
- Dispose of all opened instruments whether used or unused. This device is packaged and sterilized for single use only.
- Reuse and improper reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness, or death.
- Reuse and improper reprocessing or resterilization of single-use devices may create a risk of contamination and/or cause infection or cross-infection, including, but not limited to, the transmission of infectious diseases. Contamination may lead to injury, illness, or death.

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Sterilization

Harmonic instruments for use with Generator G11, are supplied sterile and hermetically packaged. The instruments are sterilized by gas. The sterilization process is validated conforming the requirements of EN ISO 11135-1 (for ETO sterilization of Harmonic instruments) with SAL 10^{-6} confirmed.

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


Interpretation of symbols





	For use with Gray Hand Piece only
	Use Until Date
	Lot
	Stock number
	Sterilized by ETO
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Open here
	Manufacturer/Date of Manufacture
	Authorized Representative in the European Community
	Single Patient Use
	Authorized Representative in the USA
	Temperature
	Relative Humidity
	See Instructions For Use. (Refer to blue symbol on outer packaging.)

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Do not re-sterilize

Do not use if the package is damaged

Attention! For more information see IFU

Keep from moisture, store in a dry place

Avoid of exposure direct sunlight. Keep away from sources of heat

Compliance with requirements of European Directive 93/42/EEC

Product passed procedure of declaration of conformity

Unit Quantity

Indicates that the EEP does not contain toxic and hazardous substances or elements above the maximum concentration values, and that it is an environmentally-friendly product which can be recycled (People's Republic of China symbol for Restriction of Hazardous Substances (RoHS) in electrical and electronic products).

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Keep away from moisture, store in a dry place.

Keep out of sunlight. Keep away from heat.

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
Ultrasonic instrument for generator G11
Harmonic Focus Shears + Adaptive Tissue Technology for Open Surgery

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

Important: This package insert is designed to provide instructions for use of the Harmonic FocusShears + Adaptive Tissue Technology. It is not a reference to ligation techniques.

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a Johnson & Johnson company

Instructions for use

This IFU applies to ultrasonic Instruments for the generator GEN11, with the following modifications: :

Harmonic Focus+ Shears with Adaptive Tissue Technology for Open Surgery	HAR9F
---	-------

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Indications

The Instruments are indicated for soft tissue incisions when bleeding control and minimal thermal injury are desired.

The instrument can be used as an adjunct to or substitute for electrosurgery, lasers, and steel scalpels in general, otorhinolaryngologic (ENT), plastic, pediatric, gynecologic, urologic, exposure to orthopedic structures (such as spine and joint space) and other open procedures.

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- Significant Adverse Events
- Burn and injury of Surgeon and patient
- Tissue damage, hemorrhage, leakage

If an Adverse Event or Malfunction happens, appropriate technique or action should be used.

Device Description

The Harmonic Focus+ Shears is a sterile, single patient use instrument consisting of a soft grip scissor handle housing assembly with two hand controls (MIN for minimum power level and MAX for maximum power level). The instrument has a curved blade and clamp arm with teflon pad. The instrument is 9 cm in length with a 16 mm active blade length. The Harmonic Focus+ Shears instrument allows for the cutting and coagulation of vessels up to and including 5 mm in diameter.

Each Harmonic Focus+ Shears instrument is packaged with one sterile, single patient use, disposable Gray Torque Wrench. Use only the Gray Torque Wrench with the Harmonic Focus + Shears instrument.

The torque wrench should not be discarded until the completion of the surgical case. Do not attempt to sterilize the disposable torque wrench.

The two dashes on the instrument are intended to represent relative vessel size. The MAX button is typically used for smaller vessels where cutting speed is fastest. The MIN button is typically used in slightly larger vessels and has reduced cutting speed. It is indicated for vessels up to 5 mm in size. Adaptive Tissue Technology provides the ultrasonic generator with the ability to identify and monitor the instrument during use, which enables the ultrasonic generator to modulate and decrease its power output as well as provide audible feedback to the user as appropriate.

The Harmonic Focus+ Shears is designed for use exclusively with the ultrasonic Generator G11 (GEN11) software version 2013_1 or later and Harmonic Blue Hand Piece, packaged separately. Software revision can be found under "System Information" in the ultrasonic Generator "Settings" Menu. Refer to the ultrasonic Generator Operator's Manual for more information.

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Refer to the Instructions for Use of the Harmonic Blue Hand Piece (registration certificate number FSZ 2011/11397) and Test Tip (TTBLUE) (no medical device) for instructions regarding the Hand Piece.

Illustration and Nomenclature (Illustration 1)

1. Torque Wrench
2. Blade (active jaw)
3. Clamp Arm and Tissue Pad
4. Handle Housing
5. Finger Rings
6. MIN Hand Control (proximal)
7. MAX Hand Control (distal)
8. Hand Piece (not included, registration certificate number FSZ 2011/11397)

Materials in contact with a human body:

- (2) Blade – Titanium alloy Ti 6Al4V ELI;
- (3) Clamp arm – Aluminum alloy grade 7075
- (3) Tissue Pad – Polytetrafluorethylene grade 7A
- (4) Handle Housing – Polycarbonate grade Calibre 5201-8

Operation conditions for surgical field:

Temperature: from +15 to +27°C
Relative humidity:
from 30 to 75% (without condensation)

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Storage and transportation Conditions

Temperature: from -22°C to +60°C

Relative Humidity: from 10 to 80% (without condensation)

Store away from direct sunlight, heat or freezing.

Base specifications

Specifications	Value
Clamp arm grip force, N	Not less 10.9
Targeted vessel max. diameter, mm	5
Overall sizes (A x B x C x D), mm	192 x 15,2 x 94 x 28,2
Working temperature range of blades ¹ , °C ²	40 – 300
Maximum deflection angle of the clamping jaw, °	27°
Length of active jaw, mm	16
Position of working parts at their closing	The blade should make contact along the length of the tissue pad when the working parts are closed
Sizes for Torque Wrench, mm	53,3 x 76,2 x 27,9
Mass of Torque Wrench, g	27
Requirements to the working surfaces and edges of medical instrument	Instrument surfaces shall be free of scratches, nicks, or blemishes
Requirements to the tightness (impermeability) of working parts for devices	Working parts are sealed
Requirements to the stability of the working part in biologic fluids	Working parts are resistant to blood
Mass, g	54

1-The temperature varies depending upon operating conditions

2-

40 - 120 ° C – when patient tissue is between active and passive jaws

100 - 300 ° C - active and passive jaws are closed without tissue

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<i>Specifications</i>	<i>Value</i>
Class of protection against electric shock	Class I
Type of working part	Type CF
Operating mode	Continuous duty
Time of establishment of the operating mode	Max. 10 sec
Max. output power, W	35
Max. output voltage, V	150
Max. loop resistance, Ohm	375
Maximum current, mA	400
Blade frequency, Hz	55500
Max. blade amplitude, μm	73

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Characteristics of ultrasonic electrosurgical generator GEN11 for connections of instruments

Class of protection against electric shock	Class I
Type of working part	CF type
Input characteristics	100 - 240 V ~, 50/60 Hz, 500 VA.

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Instructions for use Verify compatibility of all instruments and accessories prior to using the instrument (refer to Warnings and Precautions).

The Hand Piece and Test Tip (no medical devices) packaged separately are shipped non-sterile and must be sterilized per the insert instructions prior to each use

Preventive inspection

Before use be sure to check the tool or any other used medical devices for the absence of any defects or damage, e.g., rough surfaces, nicks or blemishes.

Assembly

- 1 Using aseptic technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field
- 2 While holding the Hand Piece, attach the instrument by rotating it onto the Hand Piece in a clockwise rotation as viewed from the distal end of the instrument (finger tight only) (Illustration 2).
- 3 Use the Gray Torque Wrench to tighten the instrument onto the Hand Piece. Turn the wrench clockwise while holding the Hand Piece until

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

it clicks twice, indicating that sufficient torque has been applied to secure the instrument (Illustration 3). To ensure proper assembly, do not grip the instrument handle while applying torque with the Torque Wrench.

Caution: Do not torque the instrument by hand or damage may occur to the Hand Piece. Do not use any means other than the Torque Wrench to attach or detach the instrument from the Hand Piece.

- 4 Remove the Torque Wrench from the instrument. Do not discard the disposable Torque Wrench until the completion of the surgical case. The Torque Wrench is used for removal of the instrument from the Hand Piece following the procedure (Illustration 4). In the event the Torque Wrench falls out of the sterile field, replace with a sterile Gray Torque Wrench. Do not re-sterilize the disposable Torque Wrench.
Caution: Take care to avoid damage to the shears when removing the Torque Wrench from the instrument.

- 5 The second activation tone can be turned off under the "Settings" Menu on the ultrasonic generator G11 (GEN11). See ultrasonic Generator Operator's Manual for more information.
This will deactivate the second activation tone only; this will not affect the Adaptive Tissue Technology's modulation and decrease of power output.

Operation

Refer to a compatible ultrasonic Generator User Manual for hand piece attachment and system operation instructions.

- 1 Connect the assembled Hand Piece and instrument to the ultrasonic generator and turn the ultrasonic generator power on. Do not turn the ultrasonic generator power on before the Hand Piece and instrument are connected to the ultrasonic generator.

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- 2 Select the desired variable or minimum power level using the INCREASE and DECREASE buttons on the ultrasonic generator.
 - Minimum starting power level defaults to power level 3 (Illustration 5). For greater tissue cutting speed use a higher ultrasonic generator power level, and for greater coagulation use a lower ultrasonic generator power level. The amount of energy delivered to the tissue and resultant tissue effects are a function of many factors, including the power level selected, blade characteristics, grip force, tissue tension, tissue type, pathology, and surgical technique.
 - MAX power is set at power level 5 and cannot be adjusted.
 - 3 The Harmonic Focus+ Shears instrument may be operated with either the foot switch or hand control. For foot switch or hand control function, refer to a compatible, ultrasonic Generator User Manual for further detail and setup and operation instructions.
 - 4 For optimal performance, clean the instrument blade and clamp arm throughout the procedure by activating the instrument tip in sterile saline (Illustration 6). The instrument can be wiped with a sterile moist gauze sponge to remove tissue, if necessary.
- WARNING: Do not touch the instrument to metal while activated (Illustration 7). See Warnings and Precautions.
- 5 If tissue is still visible in the clamp arm, use hemostats to remove residue (Illustration 8).

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Disassembly

- 1 Turn the ultrasonic generator OFF at the power switch.
- 2 Close the clamp arm and place the Gray Torque Wrench over the distal end of the instrument.
- 3 While holding the Hand Piece, loosen the instrument by turning the Torque Wrench counterclockwise. Continue to loosen by turning the instrument manually to completely unscrew it from the Hand piece.
- 4 Remove the Torque Wrench from the instrument. Dispose of the instrument and the Torque Wrench in an appropriate container.

Warnings and Precautions

- ATTENTION! Product modifications are not allowed! (Applicable for all versions of ultrasonic instruments for the electrosurgical ultrasound generator GEN11)
- Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- A thorough understanding of the principles and techniques involved in laser, electrosurgical, and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

- 6 The blade is ultrasonically energized when either the foot switch pedal is depressed or one of the hand controls is depressed.

- Pressing either the left foot pedal of the foot switch or the proximal hand control (MIN) on the instrument activates the selected minimum power level.
- Pressing either the right foot pedal of the foot switch or distal hand control (MAX) on the instrument activates the maximum power level.
- The ultrasonic generator provides an audible tone to indicate when the instrument blade is active.
- The ultrasonic generator changes to a second activation tone as Adaptive Tissue Technology regulates the delivery of energy.
- Thermal influences such as fluids or minimal to no tissue in the jaws may affect the presence or timing of the tone change.
- The tone change does not provide confirmation of tissue effect. When the second tone is heard, the situation should be assessed and the intended surgical action completed, such as gradual application of tension to facilitate transection.

- The secondary activation tone change is not a substitute for surgical experience.

WARNING: Avoid accidental contact with other instruments during use. Scratches on the blade may lead to premature blade failure. See Warnings and Precautions.

- 7 Close the clamp arm and insert the instrument through the incision. Use the Harmonic Focus+ Shears for dissection, grasping, coagulation, and cutting between the blade and clamp arm. Use the top of the blade with the clamp arm open if using for backcutting. **WARNING:** Keep the clamp arm open when backcutting or while the blade is active, without tissue between the blade and tissue pad, to avoid damage to the tissue pad.

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

electrical insulation or grounding is not compromised. Do not immerse instruments in liquid unless the instruments are designed and labeled to be immersed.

- Verify compatibility with generators. Harmonic Focus®+ Shears are compatible only with Ethicon Endo-Surgery Generator G11 (GEN11) software version 2013_1 or later. Software revision can be found under "System Information" in the Generator G11 (GEN11) "Settings" Menu. Refer to the Generator G11 (GEN11) Operator's Manual for more information.
- Audible high-pitched ringing, resonating from the blade or Hand Piece, are an abnormal condition and an indicator that the blade or Hand Piece is not operating properly. The ringing may be an indicator that the Hand Piece is beyond its useful life or that the blade has not been attached properly, which may result in abnormally high shaft temperatures and user or patient injury.
- In case of system failure, ensure the availability of the appropriate back-up equipment relevant to the specific procedure.
- Blood and tissue buildup between the blade and shaft may result in abnormally high temperatures at the distal end of the shaft. To prevent burn injury, remove any visible tissue buildup at the distal end of the shaft.
- As with all energy sources (Electrosurgery, Laser, or Ultrasound), there are concerns about the carcinogenic and infectious potential of the by-products, such as tissue smoke plume and aerosols. Appropriate measures such as protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.
 - Do not attempt to bend, sharpen, or otherwise alter the shape of the blade. Doing so may cause blade failure and user or patient injury.
 - To avoid user or patient injury in the event that accidental activation occurs, the instrument blade, clamp arm, and distal end of the shaft should not be in contact with the patient, drapes, or flammable materials while not in use.
 - During and following activation in tissue, the instrument blade and clamp arm may become hot. Avoid unintended contact with tissue, drapes, surgical gowns, or other unintended sites at all times.

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

- Incidental and prolonged activation against solid surfaces, such as bone, may result in blade heating and subsequent blade failure, and should be avoided.
- Avoid contact with any and all metal or plastic instruments or objects when the instrument is activated. Contact with staples, clips, or other instruments while the instrument is activated may result in cracked or broken blades, which may be identified by generator solid tone or instrument error.
- Scratches on the blade may lead to premature blade failure.
- Care should be taken not to apply pressure between the instrument blade and tissue pad without having tissue between them. Clamping the tissue pad against the active blade without tissue on the full length of the blade will result in higher blade, clamp arm and distal shaft temperatures and can result in possible damage to the instrument. If this happens, there may be a system failure signaled by a continuous tone or alert screen when either of the foot pedals or hand control buttons is depressed.
- Keep the jaws of the device open when backcutting or while the blade is active without tissue between the blade and tissue pad to avoid damage to the tissue pad and increased blade, clamp arm, and distal shaft temperatures.
- To avoid user or patient injury, do not activate an electrosurgical device in close proximity to the Harmonic instruments. The aerosols created by the activation of the Harmonic instruments in fatty tissue are potentially flammable.
- The entire exposed blade tip and any exposed blade shaft is active and will cut/coagulate tissue when the Harmonic Focus®+ Shears blade is activated. Be careful to avoid inadvertent contact between all exposed blade surfaces and surrounding tissue when using the Harmonic Focus + instrument.
- Use only the Harmonic Foot Switch, and the Blue Hand Piece, with the Focus + instrument to ensure compatibility with the Generator.
- After removing the instrument, examine the tissue for hemostasis. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.

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Relative Humidity

See Instructions For Use. (Refer to blue symbol on outer packaging.)

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Compliance with requirements of European Directive 93/42/EES

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Unit Quantity

Dispose of properly

Indicates that the Electrical and Electronic Product (EEP) contains toxic and hazardous substances or elements over the defined maximum concentration values. The names and content of the toxic and hazardous substances or elements shall be provided

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in the Chinese product instruction manual. The number "5" in the symbol refers to the environment-friendly use period; after their environment-friendly use period expires, they shall be put into the recycling system (People's Republic of China symbol for Restriction of Hazardous Substances (RoHS) in electrical and electronic products)

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-  Keep away from moisture, store in a dry place.
-  Keep out of sunlight. Keep away from heat.
-  Temperature
-  Relative Humidity

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Ultrasonic instrument for generator G11
**Harmonic Focus Long Shears + Adaptive Tissue Technology for
Open Surgery**

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Failure to properly follow the instructions may lead to serious surgical consequences, such as failure to ligate.

Important: This package insert is designed to provide instructions for use of the Harmonic Focus Long Shears + Adaptive Tissue Technology. It is not a reference to ligation techniques.

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Harmonic Focus+ Long Shears with Adaptive Tissue Technology for Open Surgery	HAR17F
--	--------

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Indications

The Harmonic Focus Long Shears + Adaptive Tissue Technology are indicated for soft tissue incisions when bleeding control and minimal thermal injury are desired. The instrument can be used as an adjunct to or substitute for electrosurgery, lasers, and steel scalpels in general, otorhinolaryngologic (ENT), plastic, pediatric, gynecologic, urologic, exposure to orthopedic structures (such as spine and joint space) and other open procedures.

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The following can be considered as the Adverse Event or Malfunction which may happen with use of this device.

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Significant Adverse Events

- Burn and injury of Surgeon and patient

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Device Description

The HARMONIC FOCUS+ Long Shears is a sterile, single patient use instrument consisting of a soft grip scissor handle housing assembly with two hand controls (MIN for minimum power level and MAX for maximum power level). The instrument has a curved blade and clamp arm with teflon pad. Measured from the blade tip to the MAX hand control power button, the instrument is 17 cm in length with a 16 mm active blade length. The Harmonic Focus+ Long Shears instrument allows for the cutting and coagulation of vessels up to and including 5 mm in diameter.

Each HARMONIC FOCUS+ Long Shears instrument is packaged with one sterile, single patient use, disposable purple torque wrench. Use only the purple torque wrench with the HARMONIC FOCUS+ Long Shears instrument. The torque wrench should not be discarded until the completion of the surgical case. Do not attempt to sterilize the disposable torque wrench.

The two dashes on the instrument are intended to represent relative vessel size. The MAX button  is typically used for smaller vessels where cutting speed is fastest. The MIN button  is typically used in slightly larger vessels and has reduced cutting speed. It is indicated for vessels up to 5 mm in size. Adaptive Tissue Technology provides the generator with the ability to identify and monitor the instrument during use, which enables the generator to modulate and decrease its power output as well as provide audible feedback to the user as appropriate.

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The Harmonic Focus+ Shears is designed for use exclusively with the Generator G11 (GEN11) software version 2013_1 or later and Harmonic Blue Hand Piece, packaged separately. Software revision can be found under "System Information" in the Generator G11 (GEN11) "Settings" Menu. Refer to the Generator G11 (GEN11) Operator's Manual for more information.

Refer to the Instructions for Use of the Harmonic Blue Hand Piece and Test Tip (TTBLUE) (no medical device) for instructions regarding the hand piece.

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1. Torque Wrench
2. Blade (active jaw)
3. Clamp Arm and Tissue Pad
4. Handle Housing
5. Finger Rings
6. Finger Ring Rests
7. MIN Hand Control (proximal)
8. MAX Hand Control (distal)
9. Hand Piece (not included, registration certificate number FSZ 2011/11397)

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- (2) Blade – Titanium alloy grade Ti 6Al4V ELI
- (3) Clamp arm – Aluminum alloy grade 7075
- (3) Tissue Pad – Polytetrafluorethylene grade 7A
- (4) Handle Housing – Polycarbonate grade Calibre 5201-8

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Temperature: from +15 to +27°C
Relative humidity:

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Basic specifications

Specifications	Value
Clamp arm grip force, N	no less 13,7
Targeted vessel max. diameter, mm	5
Overall sizes (A x B x C x D), mm	275 x 16 x 70 x 28
Working temperature range of blades ¹ , °C ²	40 – 300
Maximum deflection angle of the clamping jaw, °	27°
Length of active jaw, mm	16
Position of working parts at their closing	The blade should make contact along the length of the tissue pad when the working parts are closed
Sizes for Torque Wrench, mm	10,67 x 45,7 x 27,9
Mass of Torque Wrench, g	35
Requirements to the working surfaces and edges of medical instrument	Instrument surfaces shall be free of scratches, nicks, or blemishes
Requirements to the tightness (impermeability) of working parts for devices	Working parts are sealed

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Requirements to the stability of the working part in biologic fluids	Working parts are resistant to blood
Mass, g	70

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Class of protection against electric shock	Class I
Type of working part	CF type
Input characteristics	100 - 240 V ~, 50/60 Hz, 500 VA.
Output characteristics	150 VAC RMS maximum (unless otherwise specified in the instrument IFU)

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Instructions for Use

Verify compatibility of all instruments and accessories prior to using the instrument (refer to Warnings and Precautions).

The hand piece and test tip (no medical device), packaged separately, are shipped non-sterile and must be sterilized per the insert instructions prior to each use.

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Assembly

- 1 Using aseptic technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field.
- 2 While holding the hand piece, attach the instrument by rotating it onto the hand piece in a clockwise rotation as viewed from the distal end of the instrument (finger tight only) (Illustration 2).

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Operation

Refer to a compatible Harmonic Generator User Manual for hand piece attachment and system operation instructions.

- 1 Connect the assembled hand piece and instrument to the generator and turn the generator power on. Do not turn the generator power on before the hand piece and instrument are connected to the generator.
- 2 Select the desired variable or minimum power level using the INCREASE and DECREASE buttons on the generator.
 - Minimum starting power level defaults to power level 3 (Illustration 5). For greater tissue cutting speed use a higher generator power level, and for greater coagulation use a lower generator power level. The amount of energy delivered to the tissue and resultant tissue effects are a function of many factors, including the power level selected, blade characteristics, grip force, tissue tension, tissue type, pathology, and surgical technique.
 - MAX power is set at power level 5 and cannot be adjusted.
- 3 The HARMONIC FOCUS+ Long Shears instrument may be operated with either the foot switch or hand control. For foot switch or hand control function, refer to a compatible Harmonic Generator User Manual for further detail and setup and operation instructions.

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- 4 For optimal performance, clean the instrument blade and clamp arm throughout the procedure by activating the instrument tip in sterile saline (Illustration 6). The instrument can be wiped with a sterile moist gauze sponge to remove tissue, if necessary.

WARNING: Do not touch the instrument to metal while activated (Illustration 7). See **Warnings and Precautions**.

- 5 If tissue is still visible in the clamp arm, use hemostats to remove residue (Illustration 8).
- 6 The blade is ultrasonically energized when either the foot switch pedal is depressed or one of the hand controls is depressed.
- Pressing either the left foot pedal of the foot switch or the proximal hand control (MIN) on the instrument activates the selected minimum power level.
 - Pressing either the right foot pedal of the foot switch or distal hand control (MAX) on the instrument activates the maximum power level.
 - The generator provides an audible tone to indicate when the instrument blade is active.
 - The generator changes to a second activation tone as Adaptive Tissue Technology regulates the delivery of energy.
 - Thermal influences such as fluids or minimal to no tissue in the jaws may affect the presence or timing of the tone change.
 - The tone change does not provide confirmation of tissue effect. When the second tone is heard, the situation should be assessed and the intended surgical action completed, such as gradual application of tension to facilitate transection.
 - The secondary activation tone change is not a substitute for surgical experience. **WARNING:** Avoid accidental contact with other instruments during use. Scratches on the blade may lead to premature blade failure.

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- 7 Close the clamp arm and insert the instrument through the incision. Use the HARMONIC FOCUS+ Long Shears for dissection, grasping, coagulation, and cutting between the blade and clamp arm. Use the top of the blade with the jaws open if using backcutting.

Caution: Do not use finger ring rests to extend reach during procedure, as this could result in poor vessel sealing and unstable positioning of the instrument (Illustration 9).

WARNING: Keep the jaws of the device open when backcutting or while the blade is active without tissue between the blade and tissue pad, to avoid damage to the tissue pad, and increased blade, clamp arm, and distal shaft temperatures.

Disassembly

- 1 Turn the generator OFF at the power switch.
- 2 Close the clamp arm and place the purple torque wrench over the distal end of the instrument.
- 3 While holding the hand piece, loosen the instrument by turning the torque wrench counterclockwise. Continue to loosen by turning the instrument manually to completely unscrew it from the hand piece.

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- 4 Remove the torque wrench from the instrument. Dispose of the instrument and the torque wrench in an appropriate container.

Warnings and Precautions

- **ATTENTION!** Product modifications are not allowed! (Applicable for all versions of ultrasonic instruments for the electrosurgical ultrasound generator GEN11)
- Federal (USA) law restricts this device to sale by or on the order of a physician.
- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- A thorough understanding of the principles and techniques involved in laser, electrosurgical, and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that electrical insulation or grounding is not compromised. Do not immerse instruments in liquid unless the instruments are designed and labeled to be immersed.
- Verify compatibility with generators. HARMONIC FOCUS+ Long Shears are compatible only with Ethicon Endo-Surgery Generator G11 (GEN11) software version 2013_1 or later. Software revision can be found under "System Information" in the Generator G11 (GEN11) "Settings" Menu. Refer to the Generator G11 (GEN11) Operator's Manual for more information.
- Audible high-pitched ringing, resonating from the blade or Hand Piece, are an abnormal condition and an indicator that the blade or hand piece is not operating properly. The ringing may be an

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indicator that the hand piece is beyond its useful life or that the blade has not been attached properly, which may result in abnormally high shaft temperatures and user or patient injury.

- In case of system failure, ensure the availability of the appropriate back-up equipment relevant to the specific procedure.
- Blood and tissue buildup between the blade and shaft may result in abnormally high temperatures at the distal end of the shaft. To prevent burn injury, remove any visible tissue buildup at the distal end of the shaft.
- As with all energy sources (Electrosurgery, Laser, or Ultrasound), there are concerns about the carcinogenic and infectious potential of the by-products, such as tissue smoke plume and aerosols. Appropriate measures such as protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.
- Do not attempt to bend, sharpen, or otherwise alter the shape of the blade. Doing so may cause blade failure and user or patient injury.
- To avoid user or patient injury in the event that accidental activation occurs, the instrument blade, clamp arm, and distal end of the shaft should not be in contact with the patient, drapes, or flammable materials while not in use.
- During and following activation in tissue, the instrument blade, clamp arm, and distal portion of the shaft may become hot. Avoid unintended contact with tissue, drapes, surgical gowns, or other unintended sites at all times.
- Incidental and prolonged activation against solid surfaces, such as bone, may result in blade heating and subsequent blade failure, and should be avoided.
- Avoid contact with any and all metal or plastic instruments or objects when the instrument is activated. Contact with staples, clips, or other instruments while the instrument is activated may result in cracked or broken blades, which may be identified by generator solid tone or instrument error.
- Avoid accidental contact with other instruments during use. Scratches on the blade may lead to premature blade failure.

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- Care should be taken not to apply pressure between the instrument blade and tissue pad without having tissue between them. Clamping the tissue pad against the active blade without tissue on the full length of the blade will result in higher blade, clamp arm and distal shaft temperatures and can result in possible damage to the instrument. If this happens, there may be a system failure signaled by a continuous tone or alert screen when either of the foot pedals or hand control buttons is depressed.
- Keep the jaws of the device open when backcutting or while the blade is active without tissue between the blade and tissue pad, to avoid damage to the tissue pad and increased blade, clamp arm, and distal shaft temperatures.
- To avoid user or patient injury, do not activate an electrosurgical device in close proximity to the Harmonic® instruments. The aerosols created by the activation of the Harmonic instruments in fatty tissue are potentially flammable.
- The entire exposed blade tip and any exposed blade shaft is active and will cut/coagulate tissue when the HARMONIC FOCUS+ Long Shears blade is activated. Be careful to avoid inadvertent contact between all exposed blade surfaces and surrounding tissue when using the HARMONIC FOCUS+ Long instrument.
- Use only the Harmonic® Foot Switch, and the Blue Hand Piece, with the HARMONIC FOCUS+ Long Shears instrument to ensure compatibility with the Generator.
- After removing the instrument, examine the tissue for hemostasis. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.
- The sealing performance of this instrument has not been assessed on atherosclerotic vessels. Exercise caution when transecting these vessels as they may not seal properly.
- Minimum starting power level defaults to power level 3.
- Successful hemostasis may require adjunct measures when HARMONIC FOCUS+ Long Shears instruments are used on solid organs. Due to the difficulty of visualizing internal structures, proceed slowly and do not attempt to transect large masses of tissue in one activation. Avoid the division of large vascular/biliary

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bundles when using the HARMONIC FOCUS+ Long Shears instrument under these conditions.

- Products manufactured or distributed by companies not authorized by Ethicon Endo-Surgery may not be compatible with the HARMONIC FOCUS+ Long Shears instrument. Use of such products may lead to unanticipated results and possible injury to the user or patient.
- Do not torque the instrument by hand or damage may occur to the hand piece. Do not use any means other than the torque wrench to attach or detach the instrument from the hand piece.
- Take care to avoid damage to the shears when removing the torque wrench from the instrument.
- Take care to avoid injury from the blade tip while removing the torque wrench from the instrument.
- Do not clean the instrument with abrasives.
- Instruments or devices, which come into contact with bodily fluids, may require special disposal handling to prevent biological contamination.
- Dispose of all opened instruments whether used or unused. This device is packaged and sterilized for single use only.
- Reuse and improper reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness, or death.
- Reuse and improper reprocessing or resterilization of single-use devices may create a risk of contamination and/or cause infection or crossinfection, including, but not limited to, the transmission of infectious diseases. Contamination may lead to injury, illness, or death.

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Sterilization

Harmonic instruments for use with Generator G11, are supplied sterile and hermetically packaged. The instruments are sterilized by gas

The sterilization process is validated conforming the requirements of EN ISO 11135-1 (for ETO sterilization of Harmonic instruments) with SAL 10^{-6} confirmed.

For single use only. Do not attempt to re-sterilize and re-use the instrument.

How Supplied

The instruments are supplied sterile for single patient use. Each instrument is shipped with one sterile, single use, purple disposable torque wrench. Discard the instrument and torque wrench after use.

The box also includes IFU.

The cardboard packaging (box) contents:
Instrument – 6 p.
Purple Torque Wrench -6 p.

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elements over the defined maximum concentration values. The names and content of the toxic and hazardous substances or elements shall be provided in the Chinese product instruction manual. The number "5" in the symbol refers to the environment-friendly use period; after their environment-friendly use period expires, they shall be put into the recycling system (People's Republic of China symbol for Restriction of Hazardous Substances (RoHS) in electrical and electronic products).

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
Ultrasonic instrument for generator G11:
Harmonic Synergy Dissecting Hook with Torque Wrench
Harmonic Synergy Curved Blade with Torque Wrench

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
Important: This package insert is designed to provide Instructions For Use of the Harmonic Synergy instruments and Gray Torque Wrench. It is not a reference to surgical techniques.

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
Harmonic Synergy Dissecting Hook with Torque Wrench	SNGCB
Harmonic Synergy Curved Blade with Torque Wrench	SNGHK

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Indications

The Harmonic Synergy instruments are indicated for soft tissue incisions when bleeding control and minimal thermal injury are desired.

The Dissection Hook Harmonic Synergy instruments can be used as an adjunct to or substitute for electrosurgery, lasers, and steel scalpels in general, plastic, gynecologic, exposure to orthopedic structures (such as spine and joint space), ENT (Ears, Nose, Throat), including tissues of the soft palate, oral structures, and oropharyngeal airway, and thoracic surgery, including mobilization of the Internal Mammary Artery (IMA).

The Curved Blade Harmonic Synergy can be used as an adjunct to or substitute for electrosurgery, lasers, and steel scalpels in general, plastic, gynecologic, exposure to orthopedic structures (such as spine and joint space) and thoracic surgery.

Adverse events and complications

The following can be considered as the Adverse Events or Malfunctions which may happen with use of this device.

Significant Malfunctions

- Damaged and cracked blade
- Damaged and cracked sheath
- Activation issue

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- Burn and injury of Surgeon and patient
- Tissue damage, hemorrhage, leakage
- If these Adverse Events or Malfunctions happen, appropriate technique or action should be used.

Device Description

The Harmonic Synergy Dissecting Hook and Curved Blade with Torque Wrench are sterile, single patient use devices consisting of titanium blades with non-removable, telescoping gray sheaths. Two types of blades are offered: a dissecting hook and a curved blade. Selection of the appropriate blade type is a matter of surgeon preference based on the procedure or application. The Harmonic Synergy instruments allow for the cutting and coagulation of vessels up to and including 2 mm in diameter.

A soft-touch grip surface covers two sides of the handle housing. Control buttons are located on the handle housing and include a Friction Lock Button and two hand controls, MIN for minimum power level and MAX for maximum power level. When the Friction Lock Button is depressed and held down, the surgeon can adjust the length of the telescoping gray sheath up to 1.9 inches (4.8 cm) for better access and/or rotate the sheath and blade 360°, according to surgeon preference.

Each Harmonic Synergy instrument is packaged with one sterile, single patient use, disposable Gray Torque Wrench. Use only the Gray Torque Wrench with the Harmonic Synergy instruments. The Gray Torque Wrench should not be discarded until the completion of the surgical case. In the event the Torque Wrench falls out of the sterile field, replace with a sterile Gray Torque Wrench. Do not attempt to sterilize the disposable Gray Torque Wrench.

Harmonic Synergy instruments are designed for use exclusively with a compatible ultrasonic Generator and Harmonic Blue Hand Piece, packaged separately. Refer to a compatible ultrasonic Generator User Manual before using these instruments. Refer to the Instructions for Use of the Harmonic Blue Hand Piece and Test Tip (TTBLUE) (no medical device).

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1. Dissecting Hook
2. Curved Blade
3. Gray Torque Wrench
4. Sheath
5. Friction Lock Button
6. MAX hand control (distal)
7. MIN hand control (proximal)
8. Soft Grip
9. Handle Housing
10. Blue Hand Piece (not included)

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(4) Shaft (casing) – plastic Thermocomp PDX-R-98372 EES HC manufacturing by SABIC Innovative Plastics;
(1, 2) Blade – Ti 6Al4V ELI titanium alloy.
In accordance with ISO 5832-3, ASTM F1472 and ASTM B348, ASTM F136 manufacturing by Perryman Company;
(9) Handle Housing – polycarbonate Calibre 2061-15 manufacturing by Trinseo.
Primary packaging -
Lid: Tyvek 1073B with Perfecseal CR27, Perfecseal
Blister: C-Film 73X772, Amcor Flexibles
or
Lid: FMP-521, Amcor Flexibles
Blister: CI-100, Amcor Flexibles

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Storage Conditions

Temperature: from - 20°C to 25 °C
Relative Humidity: 10% – 60%

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Transport Conditions

Temperature: from -22°C to +60°C

Relative Humidity: from 10 to 65% (without condensation)

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Basic specifications

Characteristic	Value	
	SNGCB	SNGHK
Shaft diameter, mm	8,0	
Shaft adjusted length, cm	4,8	
Shaft length, cm	from 4,2 to до 9,0	
Shaft Rotation, °	360°	
Targeted vessel max. diameter, mm	2	
Working temperature range of blades ¹ , °C	40 - 100	
Overall sizes (A) x (B) x (C), mm	187,8 x 26,3 x 21,3	190,7 x 26,3 x 21,3
Size for Dissection Hook (D1) x (D2) x (D3), mm	9,5 x 4,7 x 1,0	-
Sharpening angle of blades, °	60°	30°
Size for Curved Blade, mm	-	20,7 x 2,9 x 2,2
Sizes for Torque Wrench, mm	50,8 x 68,6 x 25,4	
Mass of Torque Wrench, g	20,0	
Requirements to the working surfaces and	Instrument surfaces shall be free of scratches, nicks, or blemishes	

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edges of medical instrument	
Requirements to the stability of the working part in biologic fluids	Working parts are sealed
Requirements to the stability of the working part in biologic fluids	Working parts are resistant to blood
Mass, g	30 29

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Specifications	Value
Class of protection against electric shock	Class I
Type of working part	Type CF
Operating mode	Continuous duty
Time of establishment of the operating mode	Max. 10 sec
Max. output power, W	5.9
Max. output voltage, V	85
Max. loop resistance, Ohm	10
Maximum current, mA	400
Blade frequency, Hz	55500
Max. blade amplitude, μm	73

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Class of protection against electric shock	Class I
Type of working part	CF type
Input characteristics	100 - 240 V ~, 50/60 Hz, 500 VA
Output characteristics	150 VAC RMS maximum (unless otherwise specified in the instrument IFU) 35 watts continuous (unless otherwise specified in the instrument IFU) 30-80 kHz (55.5 kHz unless otherwise specified in the instrument IFU)

Instructions for Use

Verify compatibility of all instruments and accessories prior to using the instrument (refer to **Warnings and Precautions**).

The Blue Hand Piece and Test Tip (no medical device), packaged separately, are shipped non-sterile and must be sterilized per the insert instructions prior to each use.

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Attachment of Instrument to Blue Hand Piece

- Using sterile technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field.
NOTE: If the Gray Torque Wrench falls off the instrument when the instrument is dispensed from

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

the packaging, the Gray Torque Wrench must be re-seated on the instrument. Slide the Gray Torque Wrench, small end first, on to the instrument over the sheath. Using slight pressure to engage and align the Gray Torque Wrench ribs on the sheath, turn the Torque Wrench clockwise until engagement is felt and the Torque Wrench is seated on the sheath. Continue with the following steps.

- 2 With the Gray Torque Wrench seated over the blade, use fingers to grasp the instrument sheath and rotate the instrument onto the Blue Hand Piece in a clockwise rotation as viewed from the distal end of the instrument (finger tight only). (Illustration 2)
NOTE: Take care to avoid injury from the blade tip when attaching the instrument to the Blue Hand Piece.
- 3 While holding the Blue Hand Piece, turn the Gray Torque Wrench clockwise until it clicks twice indicating that sufficient torque has been applied to secure the instrument to the Hand Piece. (Illustration 3)
NOTE: Hold the Blue Hand Piece and not the instrument handle housing while applying the Gray Torque Wrench.
NOTE: To avoid damage to the Blue Hand Piece, do not use any means other than the Gray Torque Wrench to attach or detach the instrument from the Hand Piece.
- 4 Remove the Gray Torque Wrench from the instrument by sliding it straight back off the sheath and blade. (Illustration 4) **Do not discard the disposable Gray Torque Wrench** until the completion of the surgical procedure.
- 5 At the end of the procedure, use the Gray Torque Wrench to remove the instrument from the Blue Hand Piece. In the event the Gray Torque Wrench falls out of the sterile field, replace it with a sterile Gray Torque Wrench. Do not re-sterilize the disposable Gray Torque Wrench.

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NOTE: Take care to avoid damage to the instrument when removing the Torque Wrench from the instrument.

NOTE: Take care to avoid injury from the blade tip while removing the Torque Wrench from the instrument.

Operation

Refer to a compatible, ultrasonic Generator User Manual for Hand Piece attachment and system operation instructions.

- 1 Connect the assembled Harmonic Blue Hand Piece and Harmonic Synergy instrument to the ultrasonic generator, and turn the ultrasonic generator power on.
Caution! **Do not turn the ultrasonic generator power on before the Blue Hand Piece and instrument are connected to the ultrasonic generator.**
- 2 Select the desired power level using the INCREASE and DECREASE buttons on the ultrasonic generator. NOTE: Minimum starting power level defaults to power level 3. (Illustration 5) For greater tissue cutting speed use a higher ultrasonic generator power level, and for greater coagulation use a lower ultrasonic generator power level. The amount of energy delivered to the tissue and resultant tissue effects are a function of many factors, including the power level selected, blade characteristics, tissue tension, tissue type, pathology, and surgical technique.
- 3 The Harmonic Synergy instrument may be operated with either the foot switch or hand control. For foot switch or hand control function, refer to a compatible ultrasonic Generator User Manual for further detail and system operation instructions.

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When in use with ultrasonic Generator GEN04:

Remove the generator from STANDBY mode by depressing the stand-by button on the generator.

Allow the generator to complete the test cycle before continuing with the surgical procedure.

- 5 The instrument blade is ultrasonically energized when either foot switch pedal is depressed or one of the hand controls is depressed. Press either the left foot pedal of the footswitch or the proximal hand control (MIN) on the instrument to activate the selected minimum power level. Press either the right foot pedal of the footswitch or distal hand control (MAX) on the instrument to activate the maximum power level.

Caution! Scratches on the blade may lead to premature instrument failure.

Caution! Avoid accidental contact with other instruments during use.

- 6 After inserting the instrument through an incision, press the blade against tissue during ultrasonic activation to cut and/or coagulate tissue under direct visualization.

Caution! Because the gray instrument sheath may become hot, measures should be taken to provide a

barrier between the sheath and tissue not intended for coagulation.

Caution! In general, sharper ultrasonically activated edges cut faster with less hemostasis, while more blunt edges and surfaces coagulate more and cut less rapidly.

For optimal performance, clean the instrument blade throughout the procedure by activating the blade tip in sterile saline (Illustration 6).

Caution! Do not touch the blade to metal while activated (Illustration 7).

Caution! Do not clean the instrument with abrasives. The instrument can be wiped with a sterile moist gauze sponge to remove tissue, if necessary.

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Blade Disassembly

- 1 When in use with ultrasonic Generator GEN04: Turn the generator **OFF** at the power switch or enter **STANDBY** mode.
- 2 Slide the Gray Torque Wrench over the instrument until the Torque Wrench is properly seated on the sheath. Loosen the instrument from the Blue Hand Piece by turning the Gray Torque Wrench counterclockwise. Continue to loosen by turning the sheath manually to unscrew it completely from the Blue Hand Piece.
NOTE: Hold the Blue Hand Piece and not the handle housing while applying the Gray Torque Wrench.
procedure, verify compatibility prior to initiation of the procedure.
- 3 Remove the Gray Torque Wrench from the instrument. Dispose of the instrument and the Gray Torque Wrench in an appropriate container.
NOTE: Take care to avoid injury from the blade tip while positioning or removing the Gray Torque Wrench.

Warnings and Precautions

- **ATTENTION!** Product modifications are not allowed! (Applicable for all versions of ultrasonic instruments for the electrosurgical ultrasound generator GEN11)
- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

- A thorough understanding of the principles and techniques involved in laser, electrosurgical, and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that electrical insulation or grounding is not compromised. Do not immerse instruments in liquid unless the instruments are designed and labeled to be immersed.
- Verify compatibility with generators. The HARMONIC SYNERGY Dissecting Hook and Curved Blade instruments are not compatible with the ULTRACISION® Generator (GEN01/GEN32).
- The Dissecting Hook and Curved Blade instruments have been designed to meet the International safety standard EN60601-1 based on an intermittent operation of 15 second on/off intervals. For activation times of longer duration and under certain fault conditions, the blade sheath may become hot. To prevent burn injury, avoid direct tissue contact with the blade sheath or take preventive measures to protect tissue that comes in contact with the sheath.
- Audible high-pitched tones resonating from the instrument or Blue Hand Piece during activation are an abnormal condition and an indicator that the instrument or Blue Hand Piece is not operating properly. The tones may be an indicator that the Blue Hand Piece is beyond its useful life or that the instrument has not been attached properly, which may result in abnormally high sheath temperatures and user or patient injury.
- In case of system failure, ensure the availability of the appropriate back-up equipment relevant to the specific procedure.
- Blood and tissue buildup between the blade and sheath may result in abnormally high temperatures at the distal end of the sheath. To prevent burn injury, remove any visible tissue buildup at the distal end of the sheath.
- As with all energy sources (Electrosurgery, Laser, or Ultrasound), there are concerns about the carcinogenic and infectious potential of the by-products, such as tissue smoke plume and aerosols. Appropriate measures such as protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.
- Do not attempt to bend, sharpen, or otherwise alter the shape of the blade. Doing so may cause instrument failure and user or patient injury. • To avoid user or patient injury in the event that accidental activation occurs,

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the instrument blade and distal end of the sheath should not be in contact with the patient, drapes or flammable materials while not in use.

- During and following activation in tissue, the instrument blade may become hot. Avoid unintended contact with tissue, drapes, surgical gowns, or other unintended sites at all times.
- Incidental and prolonged activation against solid surfaces, such as bone, may result in blade heating and subsequent instrument failure, and should be avoided.
- Avoid contact with any and all metal or plastic instruments or objects when the instrument is activated. Contact with staples, clips, or other instruments while the instrument is activated may result in cracked or broken blades which may be identified by generator solid tone or instrument error.
- To avoid user or patient injury, do not activate an active electrosurgery device in close proximity to the HARMONIC SYNERGY Dissecting Hook and Curved Blade instruments. The aerosols created by the activation of the HARMONIC SYNERGY instruments in fatty tissue are potentially flammable.
- The entire exposed blade tip and any exposed blade shaft is active and will cut/coagulate tissue when the instrument blade is activated. Be careful to avoid inadvertent contact between all exposed blade surfaces and surrounding tissue when using the HARMONIC SYNERGY instrument.
- Use only the HARMONIC Foot Switch and Blue Hand Piece to ensure that they are compatible with the Generator.
- After removing the instrument, examine the tissue for hemostasis. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.
- Minimum starting power level defaults to power level 3.
- Use caution when using HARMONIC SYNERGY instruments on solid organs. Due to the limited ability of the instruments to occlude vascular structures of this nature, hemostasis may not be predictable and may require adjunct measures for coagulation.
- Products manufactured or distributed by companies not authorized by Ethicon Endo-Surgery may not be compatible with the HARMONIC SYNERGY instruments. Use of such products may lead to unanticipated results and possible injury to the user or patient.

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Dispose of all opened instruments whether used or unused. This device is packaged and sterilized for single use only.
- Reuse and improper reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness, or death.
- Reuse and improper reprocessing or resterilization of single-use devices may create a risk of contamination and/or cause infection or crossinfection, including, but not limited to, the transmission of infectious diseases. Contamination may lead to injury, illness, or death.
- If the Gray Torque Wrench falls off the instrument when the instrument is dispensed from the packaging, the Gray Torque Wrench must be re-seated on the instrument. Slide the Gray Torque Wrench, small end first, on to the instrument over the sheath. Using slight pressure to engage and align the Gray Torque Wrench ribs on the sheath, turn the Torque Wrench clockwise until engagement is felt and the Torque Wrench is seated on the sheath.
- Take care to avoid injury from the blade tip when attaching the instrument to the Blue Hand Piece.
- Hold the Blue Hand Piece and not the instrument handle housing while applying the Gray Torque Wrench.
- To avoid damage to the Blue Hand Piece, do not use any means other than the Gray Torque Wrench to attach or detach the instrument from the Hand Piece.
- Take care to avoid damage to the instrument when removing the Torque Wrench from the instrument.
- Take care to avoid injury from the blade tip while removing the Torque Wrench from the instrument.
- Do not turn the generator power on before the Blue Hand Piece and instrument are connected to the generator.
- Minimum starting power level defaults to power level 3. (Illustration 5) For greater tissue cutting speed use a higher generator power level, and for greater coagulation use a lower generator power level. The amount of energy delivered to the tissue and resultant tissue effects are a function of many factors, including the power level selected, blade characteristics, tissue tension, tissue type, pathology, and surgical technique.

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

- Scratches on the blade may lead to premature instrument failure.
- Avoid accidental contact with other instruments during use.
- Because the gray instrument sheath may become hot, measures should be taken to provide a barrier between the sheath and tissue not intended for coagulation.
- In general, sharper ultrasonically activated edges cut faster with less hemostasis, while more blunt edges and surfaces coagulate more and cut less rapidly.
- Do not touch the blade to metal while activated (Illustration 7).
- Do not clean the instrument with abrasives. The instrument can be wiped with a sterile moist gauze sponge to remove tissue, if necessary.
- Hold the Blue Hand Piece and not the handle housing while applying the Gray Torque Wrench.
- Take care to avoid injury from the blade tip while positioning or removing the Gray Torque Wrench.

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
Interpretation of symbols

Do not use this instrument with UltraCision® Generator (GEN01/GEN32).

For use with Blue Hand Piece only.

Use Until Date

Stock number

Lot

Sterilized by ETO

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

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	Manufacturer
	Authorized Representative in the European Community
	Do not reuse
	Authorized Representative in the USA
	Temperature range
	Upper limit of temperature
	Relative Humidity
	See Instructions For Use. (Refer to blue symbol on outer packaging.)
	Do not re-sterilize
	Do not use if the package is damaged
	Attention! For more information see IFU
	Keep dry
	Keep away from heat
	Compliance with requirements of European Directive 93/42/EEC
	Product passed procedure of declaration of conformity

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Indicates that the Electrical and Electronic Product (EEP) does not contain toxic and hazardous substances or elements above the maximum concentration values, and that it is an environment-friendly product which can be recycled (People's Republic of China symbol for Restriction of Hazardous Substances (RoHS) in electrical and electronic products).

Packaging unit

The transport individual packaging includes symbols:

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
Ultrasonic instrument for generator G11
Harmonic HD 1000i Shears with Integrated ergonomic handle for Open and Endoscopic Surgery
(shaft length 20 cm or 36 cm)

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Please read all information carefully.

Failure to properly follow the instructions may lead to serious surgical consequences.

Important: This Instruction for use is designed to provide instructions for use of the HARMONIC HD 1000i Shears with integrated ergonomic handle for Open and Endoscopic Surgery, shaft length 20 cm and 36 cm. It is not a reference to surgical techniques.

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




This IFU applies to ultrasonic Instruments for the generator GEN11, with the following modifications:

HARMONIC HD 1000i Shears with built-in ergonomic handle for Open and Endoscopic Surgery, shaft length 20 cm	HARHD20
HARMONIC HD 1000i Shears with built-in ergonomic handle for Open and Endoscopic Surgery, shaft length 36 cm	HARHD36

Ultrasonic Instruments for use with the Generator G11, will hereinafter be referred to as «Instrument/Instruments»

Indications

The Instruments are indicated for soft tissue incisions when bleeding control and minimal thermal injury are desired. The instruments can be used as an adjunct to or substitute for electro surgery, lasers and steel scalpels in general, plastic, pediatric, gynecologic, urologic, thoracic, exposure to orthopedic structures (such as spine and joint space), sealing and transection of lymphatic vessels, and other open and endoscopic procedures. The instruments allow for the coagulation of vessels up to and including 7 mm in diameter, using the Energy button with Advanced Hemostasis.

Device Description

The instrument is a sterile, single patient use instrument used for dissection, grasping, coagulation, and cutting between the blade and clamp arm. It consists of an ergonomic handle with an integrated hand piece and two energy delivery buttons:

- 1) Energy button - user can adjust power levels from 1-5.
- 2) Energy button with Advanced Hemostasis - for large vessel sealing; user cannot adjust.

The instrument is available in two shaft lengths - 20 cm and 36 cm.

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An integrated audible and tactile mechanism in the handle indicates full trigger closure. The instrument has a clamp arm and coated curved blade that are designed to work through a 5 mm trocar, through a 5 mm reducer cap in a larger diameter trocar, or through an incision without the use of a trocar. The instrument shaft can be rotated continuously to facilitate visualization and access to targeted tissue. The two dashes on the instrument are intended to represent relative vessel size. The Energy button is indicated for vessels up to 5 mm in diameter. When the Energy button is used, cutting speed is the fastest. The Energy button with Advanced Hemostasis is designed for larger vessels and is indicated for vessels up to 7 mm in diameter. When the Energy button with Advanced Hemostasis is used, cutting speed is reduced and hemostasis is maximized. The instrument utilizes Adaptive Tissue Technology. This provides the generator with the ability to identify and monitor the instrument during use, which enables the generator to modulate and adjust its power output as well as provide audible feedback to the user as appropriate.

This instrument is designed for use exclusively with the Generator G11 (GEN11) software version 2016-1 or later. Software revision can be found under "System Information" in the Generator G11 (GEN11) "Settings" menu. Refer to the Generator G11 (GEN11) Operator's Manual before using this instrument. Generator G11 (GEN11) is supplied separately.

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1. Coated Blade (active jaw)
2. Clamp Arm and Tissue Pad
3. Shaft
4. Rotation Knob
5. Energy Button
6. Energy Button with Advanced Hemostasis
7. Handle
8. Trigger

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- (3) Shaft – Stainless steel grade UNS S30400
- (1) Blade – Titanium alloy grade Ti 6Al4V ELI
- (1) Blade coating – Fluoropolymer grade Black E2619B
- (2) Tissue Pad – Polytetrafluoroethylene grade 7A

- (7) Handle Trigger - glass-filled silicone lubricated Polycarbonate with Styrene Block Copolymer overmold
- (7) Activation (energy) button - Polycarbonate with Styrene Block Copolymer overmold
- (7) Activation button - glass-filled Polycarbonate with Styrene Block Copolymer overmold
- (7) Handle Shrouds – Polycarbonate grade Calibre 5201-8
- (7) Rotation Knob - glass-filled silicone lubricated Polycarbonate with Styrene Block Copolymer overmold

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Temperature: from +15 to +27°C
Relative humidity: from 30 to 75% (without condensation)

Storage Conditions

Temperature: from - 20°C to 30 °C
Relative Humidity: 10% – 65%

Transport Conditions

Temperature: from -22°C to +60°C
Relative Humidity: from 10 to 65% (without condensation)

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Base technical characteristic:

Specifications	Value	
	HARHD20	HARHD36
Shaft diameter, mm	5,0	
Shaft Rotation, °	360°	

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Shaft length, cm	20	36
Maximum Indicated Vessel Size	7	
Working temperature range of blades ¹ , °C ²	40 – 300	
Blades' temperature while activation of the improved hemostasis function control button ¹ , °C ²	40 – 300	
Overall size (A) x (B) x (C), mm	380 x 133 x 33	539 x 126 x 33
Maximum deflection angle of the clamping jaw, °	46°	
Length of active jaw (blade), mm	18	
Position of working parts at their closing	"The blade should make contact along the length of the tissue pad when the working parts are closed"	
Requirements to the working surfaces and edges of medical instrument	Instrument surfaces shall be free of scratches, nicks, or blemishes	
Requirements to the tightness (impermeability) of working parts for devices	Working parts are sealed	
Requirements to the stability of the working part in biologic fluids	Working parts are resistant to blood	
Minimum width of instrument channel, mm	5,0	
Mass, g	490	

¹ - The temperature varies depending upon operating conditions

² - 40 - 120 °C – when patient tissue is between active and passive jaws
 100 - 300 °C - active and passive jaws are closed without tissue

Specifications	Value
Class of protection against electric shock	Class I
Type of working part	Type CF
Operating mode	Continuous duty
Time of establishment of the operating mode	Max. 10 sec

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Class of protection against electric shock	Class I
Type of working part	CF type
Input characteristics	100 - 240 V ~, 50/60 Hz, 500 VA.
Output characteristics	150 VAC RMS maximum (unless otherwise specified in the instrument IFU) 35 watts continuous (unless otherwise specified in the instrument IFU) 30-80 kHz (55.5 kHz unless otherwise specified in the instrument IFU)

Instructions for Use

Verify compatibility of all simultaneously used products prior to using this instrument (refer to Warnings and Precautions).

Warning and Caution Statements

WARNING: A Warning statement indicates an operating or maintenance procedure, practice, or condition that, if not strictly observed, could result in personal injury or loss of life.

Caution: A Caution statement alerts the user of a potentially hazardous

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situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property. It may also be used to alert against unsafe practices. This includes the special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Assembly

- 1 Using sterile technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field.
- 2 While holding the Hand Piece in a vertical orientation, attach the Hand Piece to the instrument by rotating the instrument onto the Hand Piece in a clockwise rotation as viewed from the distal end of the instrument (finger tight only).
- 3 Use the torque wrench (already mounted to the shaft) to tighten the blade onto the Hand Piece. Turn the torque wrench clockwise while holding only the gray Hand Piece until it clicks twice indicating that sufficient torque has been applied to secure the blade.
Note: Do not use any other means than the torque wrench to attach or detach the instrument from the Hand Piece.
Note: Do not torque the instrument by hand without the torque wrench or damage may occur to the Hand Piece.
Note: Hold only the gray Hand Piece and not the instrument handle while applying the torque wrench. (Illustration 2).
- 4 Close the trigger. Remove the torque wrench by sliding it off of the shaft. Do not dispose of the torque wrench until the procedure is completed. The

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

torque wrench is used to remove the instrument from the Hand Piece following the procedure. (Illustration 3) Dispose of the torque wrench only after completing the procedure.

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Operation

- 1 Connect the assembled Hand Piece and instrument to the generator and turn the generator power on.

Note: MIN is the only mode that can be adjusted.

- 2 Select the desired minimum power level using the INCREASE and DECREASE buttons on the generator touchscreen.

Note: The recommended minimum starting power level is Level 3 (Illustration 4). For greater tissue cutting speed use a higher generator power level, and for greater coagulation use a lower generator power level. The amount of energy delivered to the tissue and resultant tissue effects are a function of many factors, including the power level selected, blade characteristics, grip force, tissue tension, tissue type, pathology, and surgical technique.

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- 4 The blade is ultrasonically energized when either the foot switch pedal is depressed or one of the hand control buttons is depressed. Pressing the left foot pedal of the footswitch or the lower hand control button (MIN) on the instrument activates the selected minimum power level. Pressing the right foot pedal of the footswitch or upper hand control button (MAX) on the instrument activates the maximum power level. Pressing the Advanced Hemostasis button on the instrument activates the Advanced Hemostasis. Under default settings, Advanced Hemostasis is hand-activated only.

Note: The "MIN" button can be re-assigned to Advanced Hemostasis by selecting "Convert MIN to Advanced Hemostasis" on the GEN11 "Settings" screen. Refer to the Generator G11 Operator's Manual "Settings" section. (Illustration 8) Once "Convert MIN to Advanced Hemostasis" is selected, the Advanced Hemostasis button is disabled.

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Note: The generator provides an audible tone to indicate when the instrument blade is active. The generator changes to a second activation tone as Adaptive Tissue Technology regulates the delivery of energy.

- Thermal influences such as fluids or minimal to no tissue in the jaws may affect the presence or timing of the tone change.
- The tone change does not provide confirmation of tissue effect. When the second tone is heard, the situation should be assessed and the intended surgical action completed, such as gradual application of tension to facilitate transection.
- The secondary activation tone change is not a substitute for surgical experience.

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- Avoid accidental contact with other instruments during use.

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• Do not use any means other than the torque wrench to attach or detach the instrument from the Hand Piece.

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Note: To achieve complete sealing, the trigger should be fully closed and the vessel fully contained between clamp arm and blade of device. An audible and tactile “click” indicates full trigger closure.

To achieve full closure of the jaws of the device, squeeze the plastic trigger until you feel it stop against the plastic handle (plastic to plastic) (Illustration 9). If full trigger closure is released prior to or during activation on tissue, an audible and tactile “click” is evident. Increase grip force until full trigger closure is achieved.

Note: Keep clamp arm open when using the inside bottom of the blade for backcutting. (Illustration 10)

WARNING: Do not use Advanced Hemostasis mode for procedures where energy application is desired prior to full closure of the jaws (e.g. solid

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

organs). Energy is not delivered in Advanced Hemostasis mode until the jaws are completely closed.

WARNING: During benchtop testing of vessels >5 mm, the strongest vessel seals were achieved by allowing the Advanced Hemostasis mode to completely transect the targeted vessel.

WARNING: Prolonged usage of Advanced Hemostasis Mode may cause tissue pad damage.

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Disassembly

- 1 Turn the generator OFF at the power switch.
- 2 Close the clamp arm and slide the torque wrench over the distal end and up the shaft until the wrench

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aligns with the flats on the shaft. Hold by the Hand Piece only, not the instrument handle, and loosen the instrument by turning the wrench counterclockwise. Continue to loosen by turning the rotation

knob manually to completely unscrew the instrument. Do not untorque the instrument by hand

without the torque wrench or damage may occur to the Hand Piece.

Note: Do not use any other means than the torque wrench to detach the instrument from the Hand Piece.

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Warnings and Precautions

- **ATTENTION!** Product modifications are not allowed! (Applicable for all versions of ultrasonic instruments for the electrosurgical ultrasound generator GEN11)

- Federal (USA) law restricts this device to sale by or on the order of a physician.

- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.

- Minimally invasive instruments may vary in diameter from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.

- A thorough understanding of the principles and techniques involved in laser, electrosurgical, and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that electrical insulation or grounding is not compromised. Do not immerse instruments in liquid unless the instruments are designed and labeled to be immersed.

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

- Verify compatibility with generators. Use device only with Ethicon Endo-Surgery Generator G11 (GEN11) software version 2016-1 or later. Software revision can be found under "System Information" in the Generator G11 (GEN11) "Settings" menu. Refer to the Generator G11 (GEN11) Operator's Manual for more information.
- Avoid contact with any and all metal or plastic instruments or objects when the instrument is activated. Contact with staples, clips or other instruments while the instrument is activated may result in scratches on the blade, cracked or broken blades and premature blade failure.
- Do not use the Energy button with Advanced Hemostasis when energy application is desired prior to full closure of the jaws. Energy is not delivered using the Energy button with Advanced Hemostasis until the jaws are completely closed. Using the Energy button with Advanced Hemostasis without full trigger closure may result in lack of hemostasis.
- If activation is unintentionally stopped while sealing, maintain jaw closure and reactivate. Releasing the trigger while sealing may result in lack of hemostasis.
- Keep the clamp arm open when backcutting or while the blade is active without tissue between the blade and tissue pad to avoid damage to the tissue pad and increased blade, clamp arm and distal shaft temperatures.
- During benchtop testing of vessels >5 mm, the strongest vessel seals were achieved by allowing the Advanced Hemostasis mode to completely transect the targeted vessel.
- Do not clean the blade tip with abrasives. It can be wiped with a moist gauze sponge to remove tissue, if necessary. If tissue is still visible in the clamp arm, use hemostats to remove residue, taking care not to actuate the instrument. Do not touch the instrument to the hemostats while activated. Scratches on the blade may lead to cracked or broken blades and premature blade failure.
- In case of system failure, ensure the availability of the appropriate back up equipment relevant to the specific procedure.
- Audible high-pitched ringing, resonating from the blade, is an abnormal condition and an indicator that the blade is not operating properly. This may result in abnormally high shaft temperatures and user or patient injury.
- The instruments allow for the coagulation of vessels up to and including 7 mm in diameter, using the Energy button with Advanced Hemostasis. Do not attempt to seal vessels in excess of 7 mm in diameter.

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- Blood and tissue buildup between the blade and shaft may result in abnormally high temperatures at the distal end of the shaft. To prevent burn injury, remove any visible tissue buildup at the distal end of the shaft.
- As with all energy sources (Electrosurgery, Laser, or Ultrasound), there are concerns about the carcinogenic and infectious potential of the by-products, such as tissue smoke plume and aerosols. Appropriate measures such as protective eye wear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.
- Do not attempt to bend, sharpen, or otherwise alter the shape of the blade. Doing so may cause blade failure and user or patient injury.
- To avoid user or patient injury in the event that accidental activation occurs, the instrument blade, clamp arm, and distal end of the shaft should not be in contact with the patient, drapes, or flammable materials while not in use.
- During and following activation in tissue, the instrument blade, clamp arm, and distal 7 cm of the shaft may be hot. Avoid unintended contact with tissue, drapes, surgical gowns, at all times.
- Do not introduce or withdraw the instrument with the jaws open through a trocar sleeve as this may damage the instrument.
- Care should be taken not to apply pressure between the instrument blade and tissue pad without having tissue between them. Clamping the tissue pad against the active blade without tissue on the full length of the blade will result in higher blade, clamp arm and distal shaft temperatures and can result in possible damage to the instrument. If this occurs, there may be an instrument failure, and the generator touchscreen displays a troubleshooting message.
- To avoid user or patient injury, do not activate an electrosurgical device in close proximity to the HARMONIC instruments. The aerosols created by the activation of the HARMONIC instruments in fatty tissue are potentially flammable.
- The entire exposed blade tip and any exposed blade shaft are active and will cut/coagulate tissue when the instrument blade is activated. Be careful to avoid inadvertent contact between all exposed blade surfaces and surrounding tissue when using the instrument.
- Use only the appropriate foot switch, instruments, and power cord to ensure that they are compatible with the generator.
- After removing the instrument, examine the tissue for hemostasis. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

- Successful hemostasis may require adjunct measures when HARMONIC instruments are used to transect solid organs. Due to the difficulty of visualizing internal structures, proceed slowly and do not attempt to transect large masses of tissue in one activation. Avoid the division of large vascular/biliary bundles when using the instrument under these conditions.
- Products manufactured or distributed by companies not authorized by Ethicon Endo-Surgery may not be compatible with the HARMONIC system. Use of such products may lead to unanticipated results and possible injury to the user or patient.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Incidental and prolonged activation against solid surfaces, such as bone, may result in blade heating and subsequent blade failure, and should be avoided.
- Dispose of all opened instruments whether used or unused. This device is packaged and sterilized for single use only.
- Reuse and improper reprocessing or resterilization of single use devices may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness, or death.
- Reuse and improper reprocessing or resterilization of single use devices may create a risk of contamination and/or cause infection or cross-infection, including, but not limited to, the transmission of infectious diseases. Contamination may lead to injury, illness, or death.
- Tissue pad damage may occur if device is activated without tissue in the closed jaws. Activation without tissue between the jaws will cause tissue pad degradation.
- Avoid contact with any and all metal or plastic instruments or objects such as graspers, uterine manipulators.
- Prolonged blade activation with the clamp arm closed (with or without tissue between the blade and tissue pad) may cause tissue pad damage.

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Quantity

Humidity Limitation

Refer to the instructions manual/booklet

Consult Instructions for Use

Do not re-sterilize

Do not use if package is damaged.

Indicates that the Electrical and Electronic Product (EEP) contains toxic and hazardous substances or elements over the defined maximum concentration values. The names and content of the toxic and hazardous substances or elements shall be provided in the product instruction manual. The number "5" in the symbol refers to the environment friendly use period; after their environment-friendly use period expires, they shall be put into the recycling system (People's Republic of China symbol for Restriction of Hazardous Substances (RoHS) in electrical and electronic products).

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Important: This package insert is designed to provide instructions for use of the HARMONIC 5 mm

Instrument. It is not a reference to surgical techniques.

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




HARMONIC Dissecting Hook, 5 mm diameter, 32cm length	HDH05
--	-------

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Indications

HARMONIC Dissecting Hook, 5 mm diameter, 32cm length is indicated for soft tissue incisions when bleeding control and minimal thermal injury are desired. The instrument can be used as an adjunct to or substitute for electrosurgery, lasers, and steel scalpels in general, plastic, gynecologic, exposure to orthopedic structures (such as spine and joint space), and thoracic surgery, including mobilization of the Internal Mammary Artery (IMA).

Device Description

The HARMONIC 5 mm Instrument is a sterile, single patient use instrument consisting of a titanium blade with a non-removable sheath. The Instrument blade is a dissecting hook. The working length of the dissecting hook is 320 mm.

The 5 mm Instrument must be used with the 5 mm Adaptor or Hand-Switching Adaptor and connected to the HARMONIC Hand Piece and a compatible Generator prior to use. The HARMONIC 5 mm Instruments allow for the cutting and coaptive coagulation of vessels up to and including 2 mm in diameter.

The HARMONIC 5 mm Instrument is designed for use exclusively with a compatible HARMONIC Generator and Hand Piece. Refer to a compatible Generator Manual before using this instrument.

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Atmospheric range: 700 hPa to 1060 hPa

Basic technical specifications

Specifications	Value
Shaft diameter, mm	5,0 (max value)
Shaft Rotation, °	360 degrees
Blade working length, cm	32
Shaft length, cm	32
Targeted vessel max. diameter, mm	2
Working temperature range of blades ¹ , °C	10 – 50
Overall size, mm	325 x 9 x 9
Maximum torque required for coupling with laparoscopic handle, N*m	1,36
Requirements to the working surfaces and edges of medical instrument	Instrument surfaces shall be free of scratches, nicks, or blemishes
Requirements to the tightness (impermeability) of working parts for devices	Working parts are sealed
Requirements to the stability of the working part in biologic fluids	Working parts are resistant to blood
Minimum width of instrument channel, mm	5,0
Mass, g	22

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Specifications	Value
Class of protection against electric shock	Class I
Type of working part	Type CF
Operating mode	Continuous duty
Time of establishment of the operating mode	Max. 10 sec
Max. output power, W	35
Max. output voltage, V	150
Max. loop resistance, Ohm	375
Maximum current, mA	40
Blade frequency, Hz	55500
Max. blade amplitude, μm	85

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The hook blade - Titanium alloy grade Ti 6Al4V ELI, composition:

In accordance with ISO 5832-3, ASTM F1472 and ASTM B348, ASTM F136, Dynamet Inc

Body - AISI 304 stainless steel, composition

In accordance with ASTM A959, no specific raw material supplier specified.

Body coating - polydimethylsiloxane Silastic Q7-4850, natural color manufactured by Dow Corning.

Silicone primer 3-6060 Prime Coat, natural color manufactured by Dow Corning

Insulated contact - AISI 302 stainless steel composition with accordance ASTM A313

Polycarbonate Colorcomp Eastalloy DA003 EES BK8-065 779-000 black color BK8-065, manufactured by SABIC Innovative Plastics

Polyurethane ISSC 6400-5829, black with a gloss of 8-10% produced by Integrated Surgical Sciences

Vented sleeve - Polydimethylsilox Baysilone LSR 4070 natural color manufactured by Bayer

Fender - polydimethylsiloxane Silastic Q7-4765 translucent gray manufactured by Dow Corning

Tube - Teflon 100J natural color manufactured by TFX Medical

Protective shell - Polydimethylsiloxane MED-4050 Natural color manufactured by NuSil

Primary packaging -

Lid: PTH-025 Coated 1073B Tyvek, Amcor Flexibles

Blister: Clear Eastar PETG 6763 with Dow 365 Silicone, Eastman Chemical Products

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The hand piece, adaptor, and blade wrench are reusable and shipped non-sterile. These components must be sterilized per their insert instructions prior to each use.

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a Johnson & Johnson company

Attachment

- 1 Attach the 5 mm adaptor or hand switching adaptor to the hand piece. (Illustration 2)
- 2 Using sterile technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field.
- 3 Remove the protective cap from the blade tip.
Note: Take care to avoid injury from the blade tip when removing the protective cap.
- 4 Attach the blade manually to the hand piece by turning it clockwise (finger tight only).
Note: Take care to avoid injury from the blade tip when attaching the blade to the hand piece and while sliding the blade wrench onto or off of the blade. (Illustration 3)
- 5 Use the blade wrench to tighten the blade. Slide the wrench over the blade until it stops at the Hand Piece adaptor. Turn the wrench and continue to push it gently until the wrench nose is adjacent to the adaptor. Turn the wrench clockwise until it “snaps,” indicating that sufficient torque has been applied to secure the blade. (Illustration 4)
Note: Inspect the blade wrench hub for cracks or wear before use. If damage is seen, replace the blade wrench. Before use after autoclaving, cool the blade wrench at room temperature for at least 45 minutes or soak it in room temperature sterile water for 5 minutes.
- 6 Remove the blade wrench by sliding it straight back over the blade. (Illustration 5)
Do not dispose of the reusable blade wrench. It is used for removal of the blade following the procedure.

Operation

Refer to a compatible HARMONIC Generator Manual and hand piece package insert for hand piece attachment and system operation

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

- 1 Connect the assembled hand piece, adaptor, and blade to the generator and turn the generator power on.
Note: Do not turn the generator power on before the hand piece, adaptor, and blade are connected to the generator.
- 2 Select the desired power level using the INCREASE and DECREASE buttons on the generator.
- 3 The blade is ultrasonically energized when either foot switch pedal is depressed.
Note: Scratches on the blade may lead to premature blade failure.
- Avoid accidental contact with other instruments during use.
- Do not use any other means than the blade wrench to attach or detach the blade.
- 4 After inserting the blade through a trocar or an incision, press the blade against tissue during ultrasonic activation to cut and/or coagulate tissue under direct visualization.
Note: In general, sharper ultrasonically activated edges cut faster with less hemostasis, while more blunt edges and surfaces coagulate more and cut less rapidly. The amount of energy delivered to the tissue and resultant tissue effects are a function of many factors, including the power level selected, blade characteristics, tissue tension, tissue type, pathology, and surgical technique.

Disassembly

- 1 Turn the Generator OFF at the power switch or enter STANDBY mode.
- 2 Slide the wrench over the blade to the base of the hand piece. Align the flats on the wrench with the flats on the blade. Loosen the blade by turning the wrench counterclockwise. Continue to loosen by turning the wrench or blade manually to unscrew it completely. (Illustration 6)
- 3 Remove the blade wrench by pulling it straight back over the blade. Save the blade wrench for future use.

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

Note: Take care to avoid injury from the blade tip while sliding the blade wrench onto or off of the blade.

- 4 Remove the blade and dispose of it in an appropriate container. (Illustration 7)
- 5 Remove the adaptor and save it for future use.

Warnings and Precautions

- **ATTENTION!** Product modifications are not allowed! (Applicable for all versions of ultrasonic instruments for the electrosurgical ultrasound generator GEN11)
- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- A thorough understanding of the principles and techniques involved in laser, electrosurgical, and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that electrical insulation or grounding is not compromised. Do not immerse instruments in liquid unless the instruments are designed and labeled to be immersed.
- Audible high-pitched tones, resonating from the blade or hand piece, are an abnormal condition and an indicator that the blade or hand piece is not operating properly. The tones may be an indicator that the hand piece is beyond its useful life or that the blade has not been attached properly, which may result in abnormally high sheath temperatures and user or patient injury.
- Do not use the HARMONIC blade without the proper adaptor. Failure to use the proper adaptor as described in the device description may result in user or patient burn injury.
- The blade has been designed to meet the international safety standard EN60601-1 based on an **intermittent operation of 15 second on/off intervals**. For activation time of longer duration and under certain fault conditions, the blade sheath may become hot. To prevent burn injury, avoid

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direct tissue contact with the blade sheath or take preventative measures to protect tissue that comes in contact with the sheath.

- Blood and tissue buildup between the blade and sheath may result in abnormally high temperatures at the distal end of the sheath. To prevent burn injury, remove any visible tissue buildup at the distal end of the sheath.
- As with all energy sources (Electrosurgery, Laser, or Ultrasound), there are concerns about the carcinogenic and infectious potential of the by-products, such as tissue smoke plume and aerosols. Appropriate measures such as protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.
- In case of system failure, ensure the availability of the appropriate backup equipment relevant to the specific procedure.
- Do not attempt to bend, sharpen or otherwise alter the shape of the blade. Doing so may cause blade failure and user or patient injury.
- To avoid user or patient injury in the event that accidental activation occurs, the instrument blade should not be in contact with the patient, drapes or flammable materials while not in use.
- During and following activation in tissue, the instrument blade may become hot. Avoid unintended blade contact with tissue, drapes, surgical gowns, or other unintended sites after activation.
- To avoid user or patient injury, do not activate an electrosurgical device in close proximity to the HARMONIC instrument. The aerosols created by the activation of the HARMONIC instrument in fatty tissue are potentially flammable.
- The entire exposed blade tip and any exposed blade shaft is active and will cut/coagulate tissue when the instrument blade is activated. Be careful to avoid inadvertent contact between all exposed blade surfaces and surrounding tissue when using the HARMONIC instruments.
- Use only the HARMONIC Foot Switch, Hand Piece, blade accessories, and power cord to ensure that they are compatible with the Generator.
- Products manufactured or distributed by companies not authorized by Ethicon Endo-Surgery may not be compatible with the HARMONIC System. Use of such products may lead to unanticipated results and possible injury to the user or patient.
- After removing the instrument, examine the tissue for hemostasis. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Avoid contact with any and all metal or plastic instruments or objects when the instrument is activated. Contact with staples, clips, or other instruments while the instrument is activated may result in cracked or broken blades which may be identified by generator solid tone or instrument error.
- Incidental and prolonged activation against solid surfaces, such as bone, may result in blade heating and subsequent blade failure, and should be avoided.
- Dispose of all opened instruments whether used or unused.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Cleaning and resterilization of single patient use HARMONIC 5 mm Instruments can also result in abnormally high sheath temperatures and burn injury to user or patient when the blade is activated. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

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a Johnson & Johnson company

Sterilization

The Ultrasonic instruments for generator G11, various modifications, are supplied sterile and hermetically packaged. The instrument is sterilized by 25kGy γ -irradiation.

The sterilization process is validated conforming the requirements of EN ISO 11135-1 (for EO sterilization of Harmonic instruments) and the requirements of EN ISO 11137-1 (for γ -irradiation sterilization of Harmonic instruments), with SAL 10^{-6} confirmed.

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Ultrasonic instrument for generator G11

HARMONIC Disposable Hand Switching Adaptor

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Important: This package insert is designed to provide instructions for use of the Harmonic Disposable Hand Switching Adaptor. It is not a reference to surgical techniques.

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

Indication

HARMONIC Disposable Hand Switching Adaptor is designed to expand the functions of the HARMONIC Dissecting Hook, 5 mm diameter, 32cm length.

Device Description

The Harmonic Disposable Hand Switching Adaptor is a sterile, single patient use instrument designed to add hand switching capability to the Harmonic System. The hand switch includes buttons that allow the surgeon to activate maximum (MAX) or minimum (MIN) power levels. Availability of rotation allows the surgeon to adjust the location of the switch to correspond with the blade surface being used.

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1. MIN/MAX Power Switch;
2. Hand Piece Connector.

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Casing - polycarbonate Calibre 2061-15 FC850122 white colour manufactured by Trinseo;

Button - polycarbonate Calibre 2061-15 Charcoal FC780546 manufactured by Trinseo;

Label on the casing - Ethicon Ink PG PMS 307C PG/BNT/226/TP acrylic-based thermal black manufactured by ITW Transtech America Inc.;

Label on the button - Transtech America PG/2NT acrylic-based thermal black manufactured by ITW Transtech America Inc.

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ETHICON ENDO-SURGERY, LLC

a Johnson & Johnson company

Lid: Olefin polyethylene with foam coating Perfecseal CR27 Tyvek 1073B Tyvek 1073B manufactured by Perfecseal.

Bliстер: C-Film 73X772 manufactured by Amcor Flexibles

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Before using the Harmonic Disposable Hand Switching Adaptor, inspect the package for general packaging integrity. If the packaging appears to have been opened or damaged, do not use or attempt to resterilize.

The Disposable Hand Switching Adaptor is compatible with the following Harmonic System components: Generator, Hand Piece and 5 mm instrument blades. The Hand Switching Adaptor is not compatible with Harmonic Shears, 10 mm instrument blades, and Blue Hand Piece.

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Attachment

- 1 Using sterile technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field.
- 2 Attach the Hand Switching Adaptor to the Hand Piece by pushing the adaptor onto the Hand Piece until a click is heard and the adaptor is fully seated on the hand piece. (Illustration 2)
- 3 Attach the 5 mm blade following instructions in the package insert for that instrument.
- 4 Rotate the adaptor to align the MIN/MAX power switch to the desired position relative to the 5 mm instrument blade.
- 5 To activate the blade, depress the desired power level (MIN or MAX) switch. "MIN" and "MAX" positions are labeled on the switch. Refer to a compatible Harmonic Generator User Manual for a complete description of power levels.
NOTE: For hand control function, refer to a compatible Harmonic Generator User Manual for further detail and setup and operation instructions.

Disassembly

- 1 Remove the 5 mm blade following the instructions in the package insert for that instrument.
- 2 To remove the Hand Switching Adaptor from the Hand Piece, pull the adaptor firmly away from the Hand Piece.
- 3 Remove the adaptor and dispose of it in an appropriate container.
Warning: Place the Generator in Standby mode or turn off before removing or replacing an instrument, Hand Switching Adaptor or Hand Piece.

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a Johnson & Johnson company

Warnings and Precautions

- **ATTENTION!** Product modifications are not allowed! (Applicable for all versions of ultrasonic instruments for the electrosurgical ultrasound generator GEN11)
- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- A thorough understanding of the principles and techniques involved in laser, electrosurgical, and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that electrical insulation or grounding is not compromised. Do not immerse instruments in liquid unless the instruments are designed and labeled to be immersed.
- In case of system failure, ensure the availability of the appropriate back-up equipment relevant to the specific procedure.
- To avoid user or patient injury in the event that accidental activation occurs, the instrument blade and distal end of the sheath should not be in contact with the patient, drapes, or flammable materials while not in use.
- During and following activation in tissue, the instrument blade may become hot. Avoid unintended contact with tissue, drapes, surgical gowns, or other unintended sites at all times.
- Products manufactured or distributed by companies not authorized by Ethicon Endo-Surgery may not be compatible with the Harmonic System. Use of such products may lead to unanticipated results and possible injury to the user or patient.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Dispose of all opened instruments whether used or unused. This device is packaged and sterilized for single use only.
- Reuse and improper reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury or illness.

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- Reuse and improper reprocessing or resterilization of single-use devices may create a risk of contamination and/or cause infection or cross-infection, including, but not limited to, the transmission of infectious diseases. Contamination may lead to injury, illness, or death.

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How Supplied

The instruments are supplied sterile for single patient use. Discard the instrument after use.

Note: for an explanation of the symbols on the package see the end of this IFU.

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