MAQUET

INSTRUCTION FOR USE
ACROBAT SUV Vacuum Off-Pump System


Manufacturer:
Maquet Cardiovascular LLC
45 Barbour Pond Drive Wayne, NJ 07470, USA

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
stamp


2021
New Jersey, USA

**Figure 1: ACCESSRAIL Platform with
ACTIVATOR II Drive Mechanism**

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**Figure 2.
ACROBAT SUV Vacuum Stabilizer**

1. Mount Knob
2. Mount
3. Interlinking Arm
4. Tri-slot Socket
5. Malleable Stabilizer Foot
6. Vacuum Tubing with BLUE luer lock fitting

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DESCRIPTION:

The ACROBAT SUV Vacuum Off-Pump System consists of the ACCESSRAIL Platform and the ACROBAT SUV Vacuum Stabilizer. The ACCESSRAIL Platform attaches onto a reusable ACTIVATOR II Drive Mechanism. The ACCESSRAIL Platform, along with an ACTIVATOR II Drive Mechanism, is designed to create surgical access to, and direct visualization of, the thoracic cavity through a sternotomy incision. The ACCESSRAIL Platform also allows for the organization of pericardial sutures. The ACCESSRAIL Platform accommodates the stabilizer, which consists of an interlinking arm with an articulating foot at the distal end and a mount that connects to the ACCESSRAIL Platform. The stabilizer provides local immobilization of an anastomotic site.

Note: The stabilizer mount has been designed to be compatible with most adult sternotomy retractors. The mount is designed to accommodate rectangular sections ranging in thickness from 0.180" to 0.300" (4.57 to 7.62 mm) and width from 0.770" to 0.925" (19.56 to 23.50 mm). It is the responsibility of the hospital to confirm compatibility prior to beginning procedure.

INDICATIONS:

Each ACROBAT SUV Vacuum Off-Pump System is intended for use during performance of cardiac surgical procedures through a sternotomy incision approach. The ACCESSRAIL Platform in combination with an ACTIVATOR II Drive Mechanism is used to spread the sternum, providing access and direct visualization to the thoracic cavity. The ACCESSRAIL Platform also allows for the organization of pericardial sutures. The stabilizer isolates and provides local immobilization of the target vessel on the beating heart.

CONTRAINDICATIONS:

None

WARNINGS AND PRECAUTIONS:

1. Single Use Only. Do not re-sterilize. Single use medical devices can have extremely small and/or structurally complex components. It can be difficult to remove blood, tissue, body fluids and pathogens from the device. In addition, the cleaning and re-sterilization process may compromise the materials and functionality of the device.
2. Do not use if the product sterilization barrier or its packaging is compromised.
3. Do not position the stabilizer foot over a coronary artery, newly infarcted or aneurysmal heart tissue.
4. Physicians should be properly trained to perform cardiac surgical procedures with this device prior to use.
5. As with all surgical retractors, care should be taken to use only as much retraction as necessary to provide adequate access and visualization.
6. Perform the anastomosis only when the stabilizer foot is properly seated on the epicardium and adequate stabilization of the surgical site is achieved.
7. When removing the stabilizer from the epicardium, care should be taken not to disrupt the anastomotic site.
8. Take care not to adjust or move the ACCESSRAIL Platform while pericardial sutures are engaged in the platform.
9. Suture holder features are designed for "0" size sutures.
10. The ACTIVATOR II Drive Mechanism (sold separately) is supplied non-sterile and must be cleaned and sterilized prior to each use. Never use saline solution for cleaning surgical instruments.
11. The maximum vacuum setting for the ACROBAT SUV Vacuum Stabilizer is 400 mmHg.

PREPARATION FOR USE:

Thoroughly clean and sterilize the Drive Mechanism prior to each use. Follow recommended cleaning and sterilization instructions as described in the Instructions for the ACTIVATOR II Drive Mechanism.

INSTRUCTIONS FOR USE:


1. Prior to use, close the ACTIVATOR II Drive Mechanism by turning the drive handle clockwise.
2. Assemble the ACCESSRAIL Platform onto an ACTIVATOR II Drive Mechanism, matching the correct left and right platforms to the correct drive attachment side. Ensure that the drive and platforms are fully engaged.

3. Ensure that the ACCESSRAIL Platform blades are properly seated on the sternum.
4. Once the ACCESSRAIL Platform is in place, begin spreading the sternum to the desired opening by slowly turning the drive handle counter-clockwise.
5. If using sutures, slide the sutures into the suture holder slots. Engage only one suture strand per slot to ensure proper hold.
6. Release engaged suture from the platform by concurrently pulling back and up on the suture while pulling the suture through the holder slot.
7. Ensure that the Stabilizer arm moves freely prior to installing the stabilizer.
8. Assemble the stabilizer onto the ACCESSRAIL Platform by pulling back on the mount lever and hooking the mount onto the ACCESSRAIL Platform in the desired location. Lock the mount into place by sliding the mount lever 90° to the left or right.
9. Attach the BLUE luer lock fitting of ACROBAT-SUV Stabilizer Tubing Set to the BLUE luer fitting of the Vacuum Tubing on the ACROBAT-SUV Stabilizer. Attach the other end of the tubing from the Tubing Set to a fluid collection canister (not included). Make sure the Stopcock is in the closed air position. The Tubing Set with the Filter that is connected to the top of the Canister should then connect to a regulated vacuum source. Set the Vacuum Regulator to 400 mmHg. Refer to Figure 3 for proper Vacuum Tubing and system setup.
10. Gently shape malleable feet as desired to conform to the heart without bending the feet beyond 25 degrees in any direction. Place the stabilizer foot onto the heart and apply vacuum by opening the stopcock.
11. Once positioned, secure the arm by turning the mount knob clockwise.
12. To remove the stabilizer, release the vacuum by turning the stopcock to air. Loosen the stabilizer by holding the arm with one hand and turning the knob counter-clockwise with the other hand.
13. To remove stabilizer from the ACCESSRAIL Platform, slide the mount lever to the center position and pull back.
14. Unhinge the mount to release the stabilizer from the ACCESSRAIL Platform.
15. To remove the ACCESSRAIL Platform, fully close the system by turning the drive handle clockwise. Gently remove the system from the sternum.
16. To remove the ACCESSRAIL Platform from an ACTIVATOR II Drive Mechanism, turn release latches on the platform in the direction of the arrows and pull the platform away from the Drive Mechanism.

Figure 3.

ACROBAT SUV Vacuum Off-Pump System

Maximum Vacuum Setting: 400 mmHg

HOW SUPPLIED:


The ACCESSRAIL Platform and ACROBAT SUV Vacuum Off-Pump Systems are supplied **STERILE** in an unopened and undamaged package. The ACCESSRAIL Platform and ACROBAT SUV Vacuum Off-Pump System have been sterilized using gamma irradiation and are for single use only. **DO NOT RESTERILIZE. DO NOT REUSE.**

The **ACTIVATOR II Drive Mechanism** (sold separately) is a reusable instrument that is supplied non-sterile. Thoroughly clean and sterilize the Drive Mechanism prior to each use. Follow recommended cleaning and sterilization instructions as described in the Instructions for Use for the **ACTIVATOR II Drive Mechanism**.

CAUTION: Federal (U.S.A.) law restricts this device to sale, distribution and use by or on the order of a physician.

WARRANTY:

MAQUET Cardiovascular LLC (MAQUET) warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond MAQUET's control directly affect the instrument and the results obtained from its use. MAQUET's obligation under this warranty is limited to the repair or replacement of this instrument and MAQUET shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. MAQUET neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. MAQUET assumes no liability with respect to instruments reused, reprocessed or re-sterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.

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
12133

MAQUET

INSTRUCTION FOR USE ACROBAT-i Vacuum Positioner System

Manufacturer:
Maquet Cardiovascular LLC*
45 Barbour Pond Drive Wayne, NJ 07470, USA

/Signature

Huong Tran
Regulatory Affairs Specialist II
Maquet Cardiovascular LLC

stamp

James W. Monroe 6/14/2021

JAMES W. MONROE NOTARY PUBLIC OF NEW JERSEY Commission # 50142311 My Commission Expires 11/2/2025
--

2021
New Jersey, USA

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions.

Caution: Federal (USA) law restricts this device to sale by, or on the order of, a physician.

DEVICE DESCRIPTION

The ACROBAT-i Vacuum Positioner System is designed to be used in conjunction with MAQUET Cardiovascular Stabilization Systems to provide surgical access to and exposure of coronary arteries for bypass surgery. The ACROBAT-i Positioner (Positioner) is comprised of a compliant suction cup which conforms to the heart mounted to a FLEXLINK Arm (Arm) with Swivel that can be used to position the heart as desired. A locking system tightens the FLEXLINK Arm to its desired position, and a mounting system (Mount) attaches the Positioner to the ACCESSRAIL Platform, ACTIVATOR II Drive Mechanism and/or a sternal retractor. Vacuum is applied to the Suction Cup via tubing connected to a regulated vacuum source.

Note: The ACROBAT-i Positioner Mount has been designed to be compatible with a range of sizes and profiles of adult sternal retractors. It is the responsibility of the surgical practitioner to confirm that the device is attached securely to the adult sternal retractor prior to beginning the procedure.

Contents

The ACROBAT-i Vacuum Positioner System is comprised of the following (see Figure 2):

- ACROBAT-i Positioner
- ACROBAT-i Positioner Tubing Set
- ACROBAT-i Canister


Figure 1: ACCESSRAIL Platform with ACTIVATOR II Drive Mechanism (not included with this system)

1. ACTIVATOR II Drive Mechanism
2. Drive Handle
3. Platform Release Latches
4. ACCESSRAIL Platform
5. Suture Holders


Figure 2: ACROBAT-i Positioner and ACROBAT-i Canister

1. Knob with Tubing Management System
2. Mount
3. Mount Jaw
4. FLEXLINK Arm
5. Swivel
6. Grasping Handle and Suspension Assembly
7. Suction Cup (ACROBAT-i Positioner Device)
8. Vacuum Tubing with YELLOW luer lock fitting
9. Locking Lever
10. Tubing Management System on Mount
11. ACROBAT-i Canister
12. Vacuum Canister Strap
13. Volume Full Line

INDICATIONS

Each ACROBAT-i Vacuum Positioner System is intended to expose coronary arteries during the performance of cardiac surgical procedures through a sternotomy incision approach with a MAQUET Cardiovascular Stabilization System.

CONTRAINDICATIONS

None.

WARNINGS AND PRECAUTIONS

1. Single Use Only. Do not re-sterilize. Single use medical devices can have extremely small and/or structurally complex components. It can be difficult to remove blood, tissue, body fluids and pathogens from the device. In addition, the cleaning and re-sterilization process may compromise the materials and functionality of the device.
2. Do not use if the product sterilization barrier or its packaging is compromised.
3. Physicians should be properly trained to perform cardiac surgical procedures with this device prior to use.
4. Many variables, including patient anatomy, pathology, and surgical techniques, may influence procedural outcomes. Patient and procedure selection is a responsibility of the medical professional.
5. Care should be taken to ensure that mount and Locking Lever are clear of tissue when device is positioned on sternal retractor.
6. Do not engage the ACROBAT-i Positioner over a coronary artery, newly infarcted cardiac tissue or aneurismal cardiac tissue.
7. If hemodynamic instability is experienced, gently return the heart to its resting position.
8. While in operation, the ACROBAT-i Canister should be oriented and maintained in a vertical position, to ensure proper fluid collection.
9. Do not let the fluid content of the ACROBAT-i Canister exceed the full line, to prevent potential contamination of the vacuum line and the potential loss of vacuum.
10. Care should be taken to ensure that the FLEXLINK Arm is not in contact with the myocardium, to prevent unnecessary friction on the heart.
11. Care should be taken to ensure the vacuum tubings are not pinched, collapsed, or pulled during use, for proper functioning of the device.
12. The maximum vacuum setting for the ACROBAT-i Positioner is 250 mmHg.

PREPARATION FOR USE**ACTIVATOR II DRIVE MECHANISM (not included in this system)**

The ACTIVATOR II Drive Mechanism is manufactured from stainless steel. Proper care and handling of the instrument is required to ensure that the product maintains its appearance, corrosion resistance, and proper function over time.

Thoroughly clean and sterilize the Drive Mechanism prior to each use. Follow recommended cleaning and sterilization instructions as described in the Instructions for Use for the ACTIVATOR II Drive Mechanism.

INSTRUCTIONS FOR USE

1. Prepare the ACROBAT-i Canister by looping the Canister Strap and inserting the top barbed fitting of the Canister through the holes of the Canister Strap ends (as shown in Figure 3). Attach one end of the Tubing Set with the Filter to the top barbed fitting of the ACROBAT-i Canister, ensuring that it is fully seated, to create a secure holder. Hang or otherwise secure the ACROBAT-i Canister such that it is oriented and maintained in a vertical position (as shown in Figure 3).

Figure 3: ACROBAT-i Canister

Note: A hospital-supplied vacuum trap system may be substituted for the ACROBAT-i Canister.

2. Inspect the ACROBAT-i Positioner and ensure that the FLEXLINK Arm, including the Swivel and the Suspension Assembly, moves freely prior to use.
3. Attach the ACROBAT-i Positioner onto the ACCESSRAIL Platform Blades or ACTIVATOR II Drive Mechanism by pulling the Locking Lever away from the Mount body and hooking the Mount onto the retractor in the desired location while pulling back on the Mount Jaw.
4. To secure the ACROBAT-i Positioner to the retractor lock the Mount by pushing the Locking Lever towards the Mount body into the LOCKED position. Ensure the ACROBAT-i Positioner is attached securely to the retractor once the Mount is locked into place. Refer to Figure 4a and 4b for unlocked and locked positions.

Note: The Mount may be locked into place before or after positioning the FLEXLINK Arm as desired.

Note: Ensure proper clearance for Locking Lever movement to allow adequate access to the locking lever and to prevent interference with other devices.


Figure 4a: Locking Lever in UNLOCKED position


Figure 4b: Locking Lever in LOCKED position

5. Attach the **YELLOW** luer lock fitting of the ACROBAT-i Positioner Tubing Set to the **YELLOW** luer fitting of the Vacuum Tubing on the ACROBAT-i Positioner. Attach the other end of the tubing from the Tubing Set to the bottom barbed fitting of the ACROBAT-i Canister or other fluid collection canister (not included). Make sure the Stopcock is in the closed air position. The Tubing Set with the Filter that is connected to the top of the Canister should then connect from the top of the Canister to a regulated vacuum source. Set the Vacuum Regulator to 250 mmHg. Refer to Figure 5 for proper Vacuum Tubing and system setup.


Figure 5: ACROBAT-i Vacuum Positioner System
Maximum Vacuum Setting: 250 mmHg

1. Wall Source
2. Vacuum Regulator
3. Tubing Set with Filter
4. ACROBAT-i Canister
5. Tubing Set with Stopcock
6. Stopcock
7. Filter

6. Place the Suction Cup of the ACROBAT-i Positioner onto the heart and apply a vacuum by opening the Stopcock. Ensure that a tight seal is attained such that vacuum is sufficient to carry the weight of the heart. Confirm the vacuum is set at 250 mmHg.
7. Positioning of the heart should be achieved by holding the Grasping Handle at the distal segment of the FLEXLINK Arm while cradling the heart and moving it into its desired position, adjusting the Arm and Swivel position as needed. When adjusting the position of the Swivel grasp the area of the Arm adjacent to the Swivel.
Note: The positioning of the FLEXLINK Arm and Swivel may be done before or after locking the Mount into place.
8. Once positioned, ensure the Mount is locked into place and then lock the FLEXLINK Arm by turning the Knob clockwise. When turning the Knob an audible click will be heard when the peak locking force of the FLEXLINK Arm has been reached.
9. Once positioned, confirm that the ACROBAT-i Positioner Suction Cup and heart move freely with a natural movement.
10. If desired, the Vacuum Tubing on the ACROBAT-i Positioner may be secured by pressing it into the Tubing Management System grooves on the Knob or Mount.
11. To remove the ACROBAT-i Positioner, gently support the heart with one hand so it does not fall abruptly with the release of vacuum and Suction Cup disengagement. Vacuum is released by turning the Stopcock to air.
12. Loosen the FLEXLINK Arm by rotating the Knob counter-clockwise.
13. Remove the ACROBAT-i Positioner from the ACCESSRAIL Platform blade or the ACTIVATOR II Drive Mechanism by moving the Locking Lever to the UNLOCKED position, pulling the Mount Jaw back, and unhooking the Mount from the retractor.

HOW SUPPLIED

All components of the ACROBAT-i Vacuum Positioner System are supplied **STERILE** in an unopened and undamaged package. Do not use if package is opened or damaged. The ACROBAT-i Vacuum Positioner System has been sterilized using gamma irradiation and is for single use only. **DO NOT RESTERILIZE.**

DO NOT REUSE.

The ACTIVATOR II Drive Mechanism (sold separately) is a reusable instrument that is supplied non-sterile. Thoroughly clean and sterilize the Drive Mechanism prior to each use. Follow recommended cleaning and sterilization instructions as described in the Instructions for Use for the ACTIVATOR II Drive Mechanism.

WARRANTY

MAQUET Cardiovascular LLC (MAQUET) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond MAQUET's control directly affect the instrument and the results obtained from its use. MAQUET's obligation under this warranty is limited to the repair or replacement of this instrument and MAQUET shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. MAQUET neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **MAQUET assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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1133

MAQUET

INSTRUCTION FOR USE ACROBAT-i Vacuum Stabilizer System

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stamp

James W. Monroe 6/14/2021

JAMES W. MONROE
NOTARY PUBLIC OF NEW JERSEY
Commission # 60142311
My Commission Expires 11/2/2025

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DEVICE DESCRIPTION

The ACROBAT-i Vacuum Stabilizer System is designed to be used to provide surgical access to targeted vessels for bypass surgery.

The ACROBAT-i Stabilizer (Stabilizer) consists of the FLEXLINK Arm (Arm) with Swivel, an articulating and malleable Stabilizer Foot at the distal end of the Arm and a mounting system (Mount) that attaches the Stabilizer to a retractor. A locking system tightens the FLEXLINK Arm to its desired position. The Stabilizer Foot provides local immobilization of an anatomic site.

The ACROBAT-i stabilizer is a retractor based stabilizer and may be mounted on the ACCESSRAIL Platform (Standard or Deep Blades). The ACCESSRAIL Platform attaches onto a reusable ACTIVATOR II Drive Mechanism (sold separately).

Note: The ACROBAT-i Stabilizer Mount has been designed to be compatible with a range of sizes and profiles of retractors. It is the responsibility of the surgical practitioner to confirm that the device is attached securely to the retractor prior to beginning the procedure.

Contents

The ACROBAT-i Vacuum Stabilizer System is comprised of the following (see Figure 1 and 2):

- ACROBAT-i Stabilizer
- ACROBAT-i Stabilizer Tubing Set
- ACROBAT-i Canister
- ACCESSRAIL Platform (Standard Blades)


Figure 1: ACCESSRAIL Platform with ACTIVATOR II Drive Mechanism

1. ACTIVATOR II Drive Mechanism (not included in this system)
2. Drive Handle (not included in this system)
3. Platform Release Latches
4. ACCESSRAIL Platform
5. Suture Holders


Figure 2: ACROBAT-i Stabilizer and ACROBAT-i Canister

1. Knob with Tubing Management System
2. Mount
3. Mount Jaw
4. FLEXLINK Arm
5. Swivel
6. Tri-slot Socket
7. Malleable Stabilizer Foot with Mounting Pods
8. Vacuum Tubing with BLUE luer lock fitting
9. Locking Lever
10. Tubing Management System on Mount
11. ACROBAT-i Canister
12. Canister Strap
13. Volume Full Line

INDICATIONS

Each ACROBAT-i Vacuum Stabilizer System is intended for use during performance of cardiac surgical procedures. The ACROBAT-i Stabilizer isolates and provides local immobilization of the target vessel on the beating heart.

CONTRAINDICATIONS

None.

WARNINGS AND PRECAUTIONS

1. Single Use Only. Do not re-sterilize. Single use medical devices can have extremely small and/or structurally complex components. It can be difficult to remove blood, tissue, body fluids and pathogens from the device. In addition, the cleaning and re-sterilization process may compromise the materials and functionality of the device.
2. Do not use if the product sterilization barrier or its packaging is compromised.
3. Physicians should be properly trained to perform cardiac surgical procedures with this device prior to use.
4. Many variables, including patient anatomy, pathology, and surgical techniques, may influence procedural outcomes. Patient and procedure selection is a responsibility of the medical professional.
5. Care should be taken to ensure that Mount and Locking Lever are clear of tissue when device is positioned on retractor.
6. As with all surgical retractors, care should be taken to use only as much retraction as necessary to provide adequate access and visualization.
7. Do not position the ACROBAT-i Stabilizer Foot over a coronary artery, newly infarcted cardiac tissue or aneurismal cardiac tissue.
8. Perform the anastomosis only when the Stabilizer Foot is properly seated on the epicardium and adequate stabilization of the surgical site is achieved.
9. When removing the Stabilizer Foot from the epicardium, care should be taken not to disrupt the anastomotic site.
10. While in operation, the ACROBAT-i Canister should be oriented and maintained in a vertical position, to ensure proper fluid collection.
11. Do not let the fluid content of the ACROBAT-i Canister exceed the full line, to prevent potential contamination of the vacuum line and the potential loss of vacuum.
12. Care should be taken to ensure that the FLEXLINK Arm is not in contact with the myocardium, to prevent unnecessary friction on the heart.
13. Care should be taken to ensure the vacuum tubings are not pinched, collapsed, or pulled during use, for proper functioning of the device.
14. Care should be taken to not adjust or move the ACCESSRAIL Platform while pericardial sutures are engaged in the platform.
15. Suture holder features are designed for "0" size sutures.
16. The ACTIVATOR II Drive Mechanism is supplied non-sterile and must be cleaned and sterilized prior to each use. Never use saline solution for cleaning surgical instruments.
17. The maximum vacuum setting for the ACROBAT-i Stabilizer is 400 mmHg.

INSTRUCTIONS FOR USE

1. Prepare the ACROBAT-i Canister by looping the Canister Strap and inserting the top barbed fitting through the holes of the Canister Strap ends (as shown in Figure 3). Attach one end of the Tubing Set with the Filter to the top barbed fitting of the ACROBAT-i Canister, ensuring that it is fully seated, to create a secure holder. Hang or otherwise secure the ACROBAT-i Canister such that it is oriented and maintained in a vertical position (as shown in Figure 3).

Figure 3: ACROBAT-i Canister

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2. Inspect the ACROBAT-i Stabilizer and ensure that the FLEXLINK Arm, including the Swivel and the Malleable Stabilizer Foot, moves freely prior to use.
3. Attach the ACROBAT-i Stabilizer onto the ACCESSRAIL Platform blades or retractor by pulling the Locking Lever away from the Mount body and hooking the Mount onto the retractor in the desired location while pulling back on the Mount Jaw.
4. To secure the ACROBAT-i Stabilizer to the retractor lock the Mount by pushing the Locking Lever towards the Mount body into the LOCKED position. Ensure the ACROBAT-i Stabilizer is attached securely to the retractor once the Mount is locked into place. Refer to Figure 4a and 4b for unlocked and locked positions.

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5. Attach the BLUE luer lock fitting of the ACROBAT-i Stabilizer Tubing Set to the BLUE luer fitting of the Vacuum Tubing on the ACROBAT-i Stabilizer. Attach the other end of the tubing from the Tubing Set to the bottom barbed fitting of the ACROBAT-i Canister or other fluid collection canister (not included). Make sure the Stopcock is in the closed air position. The Tubing Set with the Filter that is connected to the top of the Canister should then connect to a regulated vacuum source. Set the Vacuum Regulator to 400 mmHg. Refer to Figure 5 for proper Vacuum Tubing and system setup.


Figure 5: ACROBAT-4 Vacuum Stabilizer System

Maximum Vacuum Setting: 400 mmHg

1. Wall Source

2. Vacuum Regulator

3. Tubing Set with Filter

4. ACROBAT-4 Canister

5. Tubing Set with Stopcock

6. Stopcock

7. Filter

6. Gently shape the Malleable Stabilizer Foot as desired to conform to the heart without bending the feet beyond 25 degrees in any direction.
7. Position the Malleable Stabilizer Foot at the intended location on the heart by adjusting the position of the FLEXLINK Arm and Swivel, then place the Stabilizer Foot onto the heart and apply vacuum by opening the Stopcock. When adjusting the position of the Swivel grasp the area of the Arm adjacent to the Swivel. Confirm the vacuum is set at 400 mmHg.

Note: The positioning of the FLEXLINK Arm and Swivel may be done before or after locking the Mount into place.

8. Once positioned, ensure the Mount is locked into place and then lock the FLEXLINK Arm by turning the Knob clockwise. When turning the Knob an audible click will be heard when the peak locking force of the FLEXLINK Arm has been reached.
9. If desired, the Vacuum Tubing on the ACROBAT-i Stabilizer may be secured by pressing it into the Tubing Management System grooves on the Knob or Mount.
10. To remove the ACROBAT-i Stabilizer, release the vacuum by turning the Stopcock to air and loosen the FLEXLINK Arm by rotating the Knob counter-clockwise.
11. Remove the ACROBAT-i Stabilizer from the ACCESSRAIL Platform blade or the retractor by moving the Locking Lever to the UNLOCKED position, pulling the Mount Jaw back, and unhooking the Mount from the retractor.
12. To remove the ACCESSRAIL Platform and ACTIVATOR II Drive Mechanism or other retractor, fully close the system by turning the drive handle. Gently remove the system.

HOW SUPPLIED

ACROBAT-i Vacuum Stabilizer System

All components of the ACROBAT-i Vacuum Stabilizer System are supplied **STERILE** in an unopened and undamaged package. Do not use if package is opened or damaged. The ACROBAT-i Vacuum Stabilizer System has been sterilized using gamma irradiation and is for single use only. **DO NOT RESTERILIZE. DO NOT REUSE.**

The retractor (sold separately) is a reusable instrument that is supplied non-sterile. Thoroughly clean and sterilize the retractor prior to each use. Follow recommended cleaning and sterilization instructions as described in the Instructions for Use for the retractor.

WARRANTY

MAQUET Cardiovascular LLC (MAQUET) warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond MAQUET's control directly affect the instrument and the results obtained from its use. MAQUET's obligation under this warranty is limited to the repair or replacement of this instrument and MAQUET shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. MAQUET neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. MAQUET assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.

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MAQUET is a Trademark of MAQUET GmbH & Co. KG
All rights reserved.
Copyright 2018 MAQUET Cardiovascular LLC or its affiliates. All rights reserved.