INDIVIDUAL ACKNOWLEDGMENT

State/Commonwealth of Minnesota }
County of Ramsey } ss.

On this the 27th day of August, 2021, before me,
Danielle Jean Rasmussen, the undersigned Notary Public,
personally appeared Nick Kvalheim

☒ personally known to me – OR –

☐ proved to me on the basis of satisfactory evidence

to be the person(s) whose name(s) is/are subscribed to the within instrument, and acknowledged to me that he/she/they executed the same for the purposes therein stated.

WITNESS my hand and official seal.

[Signature]

Danielle Jean Rasmussen - Notary Public
Other Required Information (Printed Name of Notary, Residence, etc.)

OPTIONAL

Although the information in this section is not required by law, it may prove valuable to persons relying on the document and could prevent fraudulent removal and reattachment of this form to another document.

Description of Any Attached Document

Title or Type of Document: Cover Letter for IFU 510663-001 Rev B, VenaSeal

Document Date: _____ Number of Pages: _____

Signer(s) Other Than Named Above: _____

Right Thumbprint of Signer
Top of thumb here

Medtronic


Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
www.medtronic.com
tel 001-763-514-4000
fax 001-763-514-4879

Cover Letter

To Whom It May Concern,

It is hereby certified that the document attached is the VenaSeal IFU (510663-001 Rev B) for the VenaSeal™ Closure System

Yours sincerely,

Signed on behalf of Medtronic, Inc.
Nick Kvalheim
Regulatory Affairs Specialist

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1 Contents

VS-403

The VenaSeal closure system is comprised of 1 vial containing 5 mL VenaSeal adhesive and the VenaSeal delivery system.

The VenaSeal delivery system consists of the following 7 delivery tools used to deliver the VenaSeal adhesive during the procedure:

- One 5 Fr (1.8 mm), 91 cm length catheter
- One 7 Fr (2.3 mm), 80 cm length introducer
- One 5 Fr (1.8 mm), 87 cm length dilator
- One 0.10 mL application dispenser gun
- Two 14 gauge dispenser tips
- Two 3 mL (3 cc) syringes
- One 0.89 mm (0.035 in) 180 cm length straight floppy tip guidewire

2 Indications for use

The VenaSeal closure system is intended for the permanent, complete, endovascular adhesive closure of the great saphenous vein (GSV) and associated varicosities in the treatment of venous reflux disease.

3 System description

	VenaSeal adhesive 5 mL of VenaSeal adhesive (specially formulated n-butyl-2-cyanoacrylate) is supplied sterile by dry heat sterilization in a vial that has a sealed top.
	VenaSeal delivery system The VenaSeal delivery system consists of 7 delivery tools. All of the tools are single use and supplied sterile (ethylene oxide sterilization).
	Catheter The catheter is 5 Fr (1.3 mm inner diameter, 1.8 mm outer diameter), with an effective length of 91 cm, laser markings at 3 cm and at 85 cm from the tip, and high echogenic visibility. The catheter is supplied sterile (ethylene oxide sterilization).
	Introducer The introducer is 7 Fr (2.0 mm inner diameter, 2.3 mm outer diameter), with an effective length of 80 cm, an echogenic tip, and circumferential markings in 10 mm intervals along its length for measuring retraction length during the VenaSeal delivery procedure. The introducer is supplied sterile (ethylene oxide sterilization).
	Dilator The dilator is 5 Fr (1.0 mm inner diameter, 1.8 mm outer diameter), with an effective length of 87 cm. The dilator is supplied sterile (ethylene oxide sterilization).
	Dispenser gun The dispenser gun consists of a pistol type, ergonomic handle with an integrated barrel and trigger. Each depression of the trigger delivers a controlled 0.10 mL (range: 0.06 mL to 0.12 mL) amount of adhesive. The dispenser gun is supplied sterile (ethylene oxide sterilization).
	Dispenser tips The dispenser tips are each comprised of a stainless steel, 1.5 mm inner diameter, 3.8 cm length hypotube with a luer lock connector. The dispenser tips are supplied sterile (ethylene oxide sterilization).
	3 mL syringes The 3 mL syringes are 3 mL, graduated Monoject™ luer lock syringes, each with a standard threaded luer lock connector. They are supplied sterile (ethylene oxide sterilization).
	Straight floppy tip guidewire The guidewire is a 0.89 mm (0.035 in) 180 cm length straight floppy tip guidewire. It is supplied sterile (ethylene oxide sterilization).

* Third-party brands are trademarks of their respective owners.

4 Contraindications

The VenaSeal closure system is not intended for use in patients with the following conditions:

- previous hypersensitivity reactions to the VenaSeal adhesive or cyanoacrylates
- acute superficial thrombophlebitis
- thrombophlebitis migrans
- acute sepsis

5 Warnings and Precautions

5.1 Warnings

- Before using, read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.
- Due to the risk of exposure to HIV or other blood borne pathogens, health care workers should always use standard blood and body fluid precautions in the care of all patients. Strictly adhere to sterile techniques during any handling of the device.
- The VenaSeal closure system should only be used by a medical provider familiar with venous anatomy and the diagnosis and treatment of conditions affecting the venous system.
- The VenaSeal system is not designed for use in the coronary and cerebral vasculature, pulmonary vasculature, or diseased and atherosclerotic arteries.
- Do not flush the catheter with saline or fluid. Flushing will make the catheter unusable with the VenaSeal adhesive.
- Manipulate the catheter only under ultrasound imaging guidance.
- Before activating the device, use ultrasound imaging to confirm that the catheter tip is positioned in the desired location.
- Inject the VenaSeal adhesive only by hand using the VenaSeal delivery system provided. Do not use a power injector.
- Do not kink the catheter. Do not use the catheter if it is kinked.
- Use the introducer in the VenaSeal delivery system for access to the peripheral vasculature. Failure to use the introducer for access may result in damage to the vessel, or failure of the VenaSeal system to perform as intended.

5.2 Precautions

- The VenaSeal system is sterile unless the package is opened or damaged. Examine the package before using the system. If the package is damaged, do not use the system.
- The VenaSeal system is supplied **sterile** for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risk of patient infections and risk of compromised device performance.
- Before using, carefully examine the VenaSeal closure system components and verify that they and their packaging have not been damaged during shipment. If the components show any sign of damage, do not use the system.
- Do not use after the use-by date.

6 Potential adverse events

The potential adverse events (e.g., complications) associated with the use of the VenaSeal system include, but are not limited to, the following list:

- Adverse reactions to a foreign body (including, but not limited to, nonspecific mild inflammation of the cutaneous and subcutaneous tissue)
- Arteriovenous fistula
- Bleeding from the access site
- Deep vein thrombosis (DVT)
- Edema in the treated leg
- Embolization, including pulmonary embolism (PE)
- Hematoma
- Hyperpigmentation
- Hypersensitivity or allergic reactions to cyanoacrylates, such as urticaria, shortness of breath, and anaphylactic shock
- Infection at the access site
- Pain
- Paresthesia
- Phlebitis
- Superficial thrombophlebitis
- Urticaria, erythema, or ulceration may occur at the injection site
- Vascular rupture and perforation
- Visible scarring

7 Storage

The VenaSeal system should be stored in its dispenser box carton at room temperature and humidity. The kit should only be opened, and the components should only be removed from the carton, upon use of the system.

8 Directions for use

8.1 Recommended accessories

The VenaSeal closure system will require the use of the following commercially available components (not supplied):

- 5 Fr (1.7 mm) micro-access kit (or equivalent)
 - 21 gauge guidewire insertion needle
 - 0.457 mm (0.018 in) micro-access kit guidewire
 - Micro-access sheath
 - Micro-access dilator
- Flushing syringe
- Sterile saline

8.2 Preparation

Prior to use, perform a baseline ultrasound to map the incompetent vein to be treated. Use sterile technique. Use local anesthetic at the access site.

1. Review the preliminary ultrasound map of the incompetent vein to be treated. Perform a final assessment using duplex ultrasound with the subject in the supine position.
2. Under ultrasound guidance, locate the great saphenous vein (GSV) and the saphenofemoral vein junction (SFJ). Trace the GSV to determine an access site at the distal end of the GSV.
3. Prepare and drape the selected access site as required per standard procedures.
4. Administer local anesthetic at the selected access site.

8.3 Access

1. Under ultrasound guidance, access the target vein with a 21 gauge guidewire insertion needle from a 5 Fr (1.7 mm) micro-access kit or equivalent using the Seldinger technique.
2. Advance a 0.457 mm (0.018 in) micro-access kit guidewire into the vein and confirm venous access by using ultrasound.
3. Remove the guidewire insertion needle and advance the micro-access sheath and dilator assembly over the micro-access kit guidewire and into the vein.
4. Remove the micro-access kit guidewire and micro-access dilator.
5. Advance the supplied 0.89 mm (0.035 in), 180 cm straight floppy tip guidewire through the micro-access sheath and into the lumen of the vein.
6. Under real-time ultrasound, position the supplied guidewire just caudal to the saphenofemoral junction (SFJ).
7. Confirm supplied guidewire positioning using ultrasound. After confirmation, remove the micro-access sheath and advance the 7 Fr (2.3 mm) introducer and dilator assembly from the VenaSeal closure system over the guidewire to the SFJ.
8. Remove the supplied guidewire and dilator.
9. Flush the introducer through the luer injection port with sterile saline using a flushing syringe. Leave the syringe in place.
10. Under real-time ultrasound guidance, position the introducer tip 5 cm (range: 4 cm to 6 cm) caudal from the SFJ.

8.4 Assembly

1. Attach a dispenser tip to the 3 mL syringe from the VenaSeal system.
2. Use the dispenser tip and 3 mL syringe combination to extract the VenaSeal adhesive from its vial into the 3 mL syringe. Purge visible air from the 3 mL syringe. See *Table 1* for the total amount of adhesive delivered as a function of planned treatment length.
3. Once the 3 mL syringe is filled with the VenaSeal adhesive, remove the dispenser tip and connect the 3 mL syringe to the 5 Fr (1.8 mm) catheter.
4. Lock the 3 mL syringe attached to the catheter into place within the dispenser gun. Use caution to avoid kinking the catheter. *Figure 1* illustrates the assembly of the VenaSeal system.

Figure 1. Illustration of the dispenser gun, 3 mL syringe, and catheter assembly

- | | |
|-----------------|----------------|
| 1 Dispenser gun | 3 3 mL syringe |
| 2 Trigger | 4 Catheter |

5. Prime the catheter by pulling the trigger on the dispenser gun to advance the VenaSeal adhesive to within 3 cm (range: 1 cm to 5 cm) of the distal catheter tip. The distal laser mark denotes 3 cm from the catheter tip. Do not advance the adhesive all the way to the catheter tip during priming.
6. Remove the saline-filled flushing syringe from the introducer.
7. Insert the primed catheter into the introducer and advance until the proximal laser mark on the catheter is at the hub of the introducer. Use caution to avoid kinking the catheter.
8. Pull the introducer back caudal another 5 cm. Advance the catheter cephalad and lock the introducer to the catheter with the spin-lock mechanism. This action will expose the catheter tip, positioning it 5 cm (range: 4 cm to 6 cm) caudal to the SFJ. Verify that the catheter tip is 5 cm from the SFJ using ultrasound.

8.5 Procedure

8.5.1 Initial delivery of the adhesive

1. Using a transverse ultrasound plane, position the ultrasound transducer just cephalad (2 cm to 3 cm) to the catheter tip and apply adequate pressure to compress the GSV near the SFJ as depicted in *Figure 2*.

Figure 2. Ultrasound transducer position for initial injection of the VenaSeal adhesive

2. While applying compression with the ultrasound transducer, deliver an aliquot of 0.10 mL (range: 0.06 mL to 0.12 mL) of the VenaSeal adhesive in the vein by depressing the trigger of the dispenser gun once and holding it down for 3 seconds to completely inject the 0.10 mL of the VenaSeal adhesive. Immediately pull back 1 cm and deliver another aliquot of 0.10 mL of the VenaSeal adhesive in the vein with an additional trigger pull, again holding down the dispenser gun trigger for 3 seconds after depressing the trigger.
3. Following this second injection, immediately (within 3 seconds) pull back the introducer and catheter assembly 3 cm, holding the transverse compression with the ultrasound transducer for a minimum of 3 minutes. During the compression time, use a free hand caudal to the transducer to add light compression as depicted in *Figure 3*.

Figure 3. Pressure application in the VenaSeal system procedure


8.5.2 Subsequent delivery of the adhesive

The remaining length of vein that requires treatment can be treated with any combination of 3 cm, 6 cm, or 9 cm segments as described in the following procedure.


1. Locate the catheter tip under ultrasound guidance.
 2. Apply compression on the vein caudal to the previous injection and cephalad to the catheter tip with the ultrasound transducer in a transverse plane.
 3. Deliver an aliquot of 0.10 mL (range: 0.06 mL to 0.12 mL) of the VenaSeal adhesive in the vein by pulling the trigger of the dispenser gun once and holding the trigger for 3 seconds.
 4. Following the injection, immediately (within 3 seconds) pull back the introducer and catheter assembly 3 cm. See *Figure 4*.
- Note:** To treat a 6 cm segment before compression, repeat steps 3 and 4 once to dispense 2 aliquots of adhesive. See *Figure 5*.
- Note:** To treat a 9 cm segment before compression, repeat steps 3 and 4 twice to dispense 3 aliquots of adhesive. See *Figure 6*.
- Note:** The treating physician should use professional judgement to determine the exact application, considering the remaining length of vein to be treated, anatomical considerations (for example, tributaries, perforators, locations of focal dilation), and the individual patient situation.
5. While maintaining transverse compression with the ultrasound transducer, use a free hand caudal to the transducer to add light compression to the entire segment treated in *Step 3* and *Step 4* for a minimum of 30 seconds. Refer to *Figure 4*, *Figure 5*, or *Figure 6* depending on the length of segment being compressed.

Figure 4. Treatment of 3 cm segment


- 1 Delivery of adhesive and pullback
- 2 Thirty seconds of compression

Figure 5. Treatment of 6 cm segment

- 1 Delivery of 2 aliquots of adhesive with pullback
- 2 Thirty seconds of compression

Figure 6. Treatment of 9 cm segment

- 1 Delivery of 3 aliquots of adhesive with pullback
- 2 Thirty seconds of compression

6. Repeat Step 1 to Step 5 to treat the entire length of the target vein segment to a point where the 5 cm laser mark of the introducer is visible outside of the access site. Injections should all be located within the target vein. The introducer should also remain within the target vein.
- Note:** Additional injections in the same location can be given during treatment at the site of tributaries or focal dilatation.
7. After the 30 seconds of compression associated with the final injection within the target vein, unlock the spin-lock mechanism of the catheter from the introducer. Maintain the introducer position within the target vessel.
8. While holding the introducer stationary, recapture the catheter by retracting it through the introducer until the catheter's proximal laser mark is visible 1 cm to 5 cm outside of the introducer hub.
9. Remove the introducer and catheter together. Apply hand pressure as long as necessary to achieve hemostasis at the access site.
10. Confirm vein closure along the treated segment by using ultrasound.

8.6 If the primed catheter needs to be withdrawn from the introducer prior to completion, take the following steps:

1. If the adhesive has not been delivered through the distal end of the catheter (Section 8.4, Step 6 to Section 8.5.1, Step 1), perform the following steps:
 - a. Ensure that no adhesive has been delivered from the catheter tip.
 - b. If no adhesive has been delivered from the catheter tip, unscrew the catheter and withdraw it.
2. If the adhesive has been injected through the distal end of the catheter (after Section 8.5.1, Step 1), perform the following steps:
 - a. Withdraw both the introducer and the catheter by following Section 8.5.2, Step 7 to Section 8.5.2, Step 9.
 - Note:** Step 2b to Step 2i are performed outside the body.
 - b. Re-expose the tip of the catheter and re-lock the catheter to the introducer with the spin-lock mechanism.
 - c. Push the white button on top of the dispenser gun and pull back the plunger rod.
 - d. Unlock and release the syringe with the attached catheter from the dispenser gun.
 - e. Pull back gently on the syringe plunger to draw the adhesive back into the catheter.
 - f. The adhesive should be drawn at least 5 cm to 10 cm from the distal laser mark on the catheter.
 - g. Remove any residual adhesive by repeatedly wiping the tip of the catheter with sterile gauze.
 - h. After confirming that the tip is cleared of adhesive, unscrew the catheter from the introducer and withdraw the catheter.
 - i. Flush the introducer sheath with sterile saline.
3. To reintroduce the catheter into the introducer, follow Section 8.3, Step 6 to Section 8.5.1, Step 1.

8.7 Total amounts of VenaSeal adhesive expected for use in vein treatments

The amount of VenaSeal adhesive delivered is related to length of the target GSV. The VenaSeal system instructions for use describe the application of 0.10 mL of the VenaSeal adhesive every 3 cm in the target vein. In addition, the first application consists of 0.20 mL. Table 1 shows a calculation of the total expected amount of adhesive delivered as a function of target treatment length. The exact application procedure and dosage is to be determined by the medical provider.

Table 1. Reference chart for total amount of adhesive delivered as a function of planned treatment length.

Planned treatment length (cm)	Initial delivery (mL)	Subsequent number of deliveries	Total amount of adhesive (mL)
10	0.2	3	0.5
15	0.2	5	0.7
20	0.2	7	0.9
25	0.2	8	1
30	0.2	10	1.2
35	0.2	12	1.4
40	0.2	13	1.5
45	0.2	15	1.7
50	0.2	17	1.9
55	0.2	18	2
60	0.2	20	2.2
65	0.2	22	2.4
70	0.2	23	2.5

9 Disposal

Dispose of the VenaSeal system in accordance with applicable laws, regulations, and hospital procedures regarding biohazardous medical devices.

10 Disclaimer of warranty

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations

shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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
STATE OF MINNESOTA

COUNTY OF HENNEPIN

I, Ashley Sophia Peterson, Notary for the State of Minnesota, hereby certify that the attached Statement Letter-Appendix for the IFU, VenaSeal is a true copy of the original filed at Medtronic, Inc.

Dated: September 26, 2022

(stamp)

Ashley Peterson
Signature of Notary Public

Title: Notary Public
My commission expires: January 31, 2027

Medtronic


Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA


www.medtronic.com


Statement Letter

To Whom It May Concern:

Medtronic, Inc., a corporation organized under the laws of the United States, having its principal place of business at 710 Medtronic Parkway, Minneapolis, MN 55432, USA, as the legal manufacturer of the 'VenaSeal Closure System', hereby confirm that the attached document – Appendix to the Instruction For Use, is a true original document.

Erin Roquemore
International Regulatory Affairs Specialist

Date

Appendix to the Instruction for use
VenaSeal Closure System

(When the device is delivered to the Russian Federation, the information contained in this appendix to the technical manual is a priority)

NAME OF MEDICAL DEVICE

VenaSeal Closure System as part of:

- 1. VenaSeal adhesive, 5 ml – one piece.
- 2. VenaSeal delivery system:
 - 2.1. Catheter 5 Fr – one piece.
 - 2.2. Introducer 7 Fr – one piece.
 - 2.3. Dilator 5 Fr – one piece.
 - 2.4. Application dispenser gun - one piece.
 - 2.5. Gauge dispenser tips – two pieces.
 - 2.6. Syringes, 3 ml – two pieces.
 - 2.7. Straight floppy tip guidewire – one piece.

INFORMATION ON RUSSIAN STICKER:

The final packaging for the Russian market is accompanied with additional sticker.
Devices supplied in packaging marked with the following possible information:

- Name of medical device in Russian language;
- Name of the manufacturer;
- Country of origin;
- Code UPN (unique product number);
- Code CFN (customer facing number);
- Date and number of Registration Certificate;
- Warning and informing symbols and inscriptions (if applicable)

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STERILIZATION

VenaSeal Closure System as part of:	
1. VenaSeal adhesive, 5 ml - one.	Sterilized using dry heat
2. VenaSeal delivery system:	Sterilized using ethylene oxide
2.1. Catheter 5 Fr - one.	
2.2. Introducer 7 Fr - one.	
2.3. Dilator 5 Fr - one.	
2.4. Application dispenser gun - one.	
2.5. Gauge dispenser tips - two.	
2.6. Syringes, 3 ml - two.	
2.7. Straight floppy tip guidewire - one.	

REQUIREMENTS FOR THE PROTECTION OF THE ENVIRONMENT IN THE APPLICATION OF MEDICAL DEVICES. PROCEDURE FOR UTILIZATION/LIQUIDATION OF MD

After use, this product may be a potential biohazard (hazard class B in accordance with SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of territories of urban and rural settlements, for water bodies, drinking water and drinking water supply, atmospheric air, soils, residential premises, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures"). Handle and dispose of all such device are accordance with accepted medical practice and applicable local, regional and national laws and regulations.

Unused products with expired shelf life are classified as Class A medical waste.

CONDITIONS OF TRANSPORTATION AND STORAGE

Temperature	+15 to +30 °C
Relative humidity	30% -70%
Pressure	700 – 1060 hPa

OPERATING CONDITIONS

Temperature	+32°C to +42°C (Body Temperature)
Relative humidity	100% non-condensing
Pressure	Not applicable, non-hermetic product

SHELF LIFE

2 years.

PROCEDURE AND CONDITIONS FOR THE DISPOSAL OR DESTRUCTION OF A MEDICAL DEVICE

After use, this product may be a potential biohazard (hazard class B in accordance with SanPiN). Handle and dispose of all such device are accordance with accepted medical practice and applicable local, regional and national laws and regulations.

THE LIST OF NATIONAL AND INTERNATIONAL REGULATORY DOCUMENTS/STANDARDS, WHICH MATCHES MEDICAL DEVICE

EN ISO 13485, EN ISO 14971, ASTM F1929-12, ASTM F88-09, ASTM D4169-14, ASTM F2096-11, EN ISO 11607-2, EN ISO 11607-1, ASTM F1980-02, EN ISO 10993-7, EN ISO 11135-1, EN 556-1, MEDDEV 2.7.1 Rev 3, EN ISO 14630, IEC 62366, BS EN 1041, EN 980, ISO 15223-1, EN

Appendix to the Instruction for use
VenaSeal Closure System

ISO 10993-1, EN ISO 10993-3, EN ISO 10993-4, EN ISO 10993-5, ISO 10993-10, EN ISO 10993-11, EN ISO 10993-12, EN ISO 10555-1, BS EN ISO 20594-1.

CUSTOMER SERVICE

Contact for all questions, including complaints of medical devices the Authorized Representative in Russia:

Medtronic, LLC,

123112, Russia, Moscow, in.ter.c. Presnensky Municipal District, Presnenskaya nab., 10, floor 9, office III, room 41.

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© 2020 Medtronic