
Ref. : 002-21

Date: 09 December 2021

APPROVED

Cristiano Rizzo / CEO / Micromed S.p.A.

« 09th » December **2021**

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ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL RESEARCH
"MICROMED", IN VERSIONS
Version 1.2

Production: MICROMED S.P.A., Via Giotto 2, Mogliano Veneto (TV) CAP 31021, Italy

2021

Micromed Group | located in ITALY · BELGIUM · FRANCE · GERMANY · USA | www.micromedgroup.com

BRAIN QUICK

MOBERG

NEUROWERK

MYOQUICK

SleepRT

Micromed S.p.A. a socio unico
Via Giotto 2
31021 Mogliano Veneto (TV)

Phone: +39.041.5937000
Fax: +39.041.5937011
micromed.it@micromedgroup.com

P.IVA / VAT No: IT 03906850262
Reg. Impr. e CF TV03906850262
REA n.307680

Capitale sociale 1.000.000,00€
Iscri. Reg. AEE IT08020000000841

I the undersigned Nicoletta Mammucari, Notary public in Mogliano Veneto, duly enrolled in the Notarial College of the District of Treviso, who knows English language, attest the authenticity of the signature of Mr. Cristiano RIZZO, born in Venice on the 2th of March 1969, resident for the office as follows, as special attorney of the Company MICROMED SPA, with registered office in Mogliano Veneto (TV) Via Giotto n.2, tax code and registration number in the Treviso - Belluno Company Register n.03906850262, right power of attorney of 24th of July 2018, duly filed in the company register, whose personal identity, qualification and powers I, the Notary, am certain.

Mogliano Veneto, Via IV Novembre n.6/1, on the tenth of December twothousandtwentyone.

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MICROMED S.P.A.

14th February 2023

Appendix to Operational Documentation for a medical device:

**ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL RESEARCH
"MICROMED"**

Version 2

Production: MICROMED S.P.A., Via Giotto 2, Mogliano Veneto (TV) CAP 31021, Italy

2023

1. Purpose and scope.

The electroencephalographic (EEG) system for neurophysiological studies " MICROMED " is designed for graphic recording of bioelectrical signals, visual display of the electrical activity of the brain in order to diagnose neurological diseases characterized by episodic changes in EEG parameters or evoked potential (hereinafter referred to as the EEG system).

The system has 9 types of executions. Each version of the system differs in the number of EEG channels in EEG blocks, in all types of versions, without exception, identical software is installed, in addition, all components (excluding EEG blocks and amplifiers for them, since the types of system versions differ in the number of channels) of types of system versions are the same.

Scope: hospitals, clinics, hospitals, medical institutions, outpatient clinics.

Potential consumers: Qualified medical personnel (nurse, paramedic, functional diagnostics doctor) only for electroencephalography and evoked potential studies.

Attention! When working with the system and its components, as well as when assembling the system, be sure to: wear disposable gloves to avoid personnel contact with the system and the transfer of infectious agents to external parts of the system. You also need to be wearing a mask. ATTENTION! GLOVES AND MASK ARE NOT INCLUDED AND ARE NOT SUPPLIED TOGETHER WITH THE PRODUCT.

The patient is also required to wear disposable gloves (it is allowed to remove gloves only when working with sensors that are attached to the finger). This requirement was introduced in connection with the control of the situation on COVID -19.

ATTENTION! GLOVES AND MASK ARE NOT INCLUDED AND ARE NOT SUPPLIED TOGETHER WITH THE PRODUCT.

The contact parts are only: electrodes, caps and sensors, tapes.

A list of referenced documents and European Standards is given in Appendix B.

Indications for use:

- sleep disorders;
- convulsive attacks;
- identified endocrine diseases;
- craniocerebral injuries (assessment of the severity and effectiveness of the restoration of brain function after an injury). With repeated studies, the EEG helps to assess the speed and completeness of the disappearance of signs of brain dysfunction. Further treatment depends on this ;
 - pathology of the vessels of the head and neck (identified by ultrasound);
 - inflammatory diseases of the nervous system and their consequences, namely:
 - ✓ bacterial meningitis or meningoencephalitis;
 - ✓ cerebral malaria;

- ✓ cerebral toxoplasmosis;
- ✓ cytomegalovirus infection;
- ✓ HIV;
- ✓ tuberculosis;
- ✓ viral encephalitis;
- ✓ subacute sclerosing panencephalitis;
- vegetovascular dystonia;
- panic attacks, diencephalic crises;
- cerebrovascular diseases: transient ischemic attacks, stroke and its consequences, anomalies in the structure of cerebral vessels;
- frequent headaches, dizziness;
- after a neurosurgical operation;
- with any brain damage that developed before or after childbirth;
- developmental and behavioral disorders in children;
- volumetric formations of the brain;
- neoplasm is suspected. By modern standards, an EEG study can be recommended as a screening study. Due to the harmlessness, relative availability and speed of the EEG, the doctor can tell whether it is worth sending the patient for an additional (more often tomographic) study or not.

- epilepsy, to establish the areas of the brain involved in the triggering of seizures, to monitor the effect of the action of drugs, to decide on the termination of drug therapy, to determine the degree of brain dysfunction during periods when there are no seizures;

- sleep disorders of a paroxysmal nature;
- behavioral disorders of a paroxysmal nature (including without impaired consciousness);
- fainting (syncope);
- to clarify the diagnosis of patients with:
 - ✓ Headaches.
 - ✓ Vertigo.
 - ✓ Unstable blood pressure, hypertension.
 - ✓ Vegetovascular dystonia.
 - ✓ Vertebral-basilar insufficiency in cervical osteochondrosis.
 - ✓ neurotic disorders;
 - ✓ Discirculatory encephalopathy;
- to study the functional state of the brain in people in whom structural research methods (for example, MRI: magnetic resonance imaging method) show that the brain structure is without pathology, but there is a violation of brain function in the clinic (for example, with metabolic encephalopathy);
 - acute disorders of cerebral circulation (early and long-term consequences of strokes);
 - perinatal pathology of the nervous system (PP CNS in infants and young children);
 - delayed psycho-speech development of unclear origin;

EEG allows you to monitor the results of treatment in patients receiving one or another type of drug or physiotherapy. It is possible to conduct a survey of patients of any age and in any functional state of the patient (including sleep).

Contraindications for use:

The procedure has practically no contraindications for use. The only thing that can limit its implementation is the presence of extensive injuries on the surface of the head, acute infectious processes or postoperative sutures that have not healed by the time of the study. In addition, a contraindication is the presence of a history of allergic reactions to the materials used in the manufacture of the product.

Attention! The study is not contraindicated for children and pregnant women, but during these periods it is performed without functional tests.

Possible Complications

There were no side effects or any complications after the EEG .

A possible complication may be an allergic reaction to the materials used in the manufacture of the product.

Potential consumers: Qualified medical personnel (nurse, paramedic, functional diagnostics doctor) only for electroencephalography and evoked potential studies.

1.1 Contents of delivery.

ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL STUDIES MICROMED , IN THE FOLLOWING VERSIONS:

1. ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL STUDIES MICROMED , PERFORMED BY: SD LTM 32 EXPRESS , CONSISTING OF:

System software PLUS EVOLUTION - BRAIN QUICK including the following software options:

Software option for registration and viewing:

- 2.1 Software option for recording and viewing EEG, 2-channel CFM version (on request)
- 2.2 Software option for recording and viewing EEG, 8-channel CFM version (on request)
- 2.3 Software option for recording and viewing EEG, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
- 2.4 Software option for recording and viewing digital video-EEG up to 64 channels, version CLINI C (on request)
- 2.5 Software option for recording and viewing digital video-EEG, version CLINI C (on request)
- 2.6. Software option for recording and viewing digital video-EEG in Full HD format (on request)
- 2.7 Software option for recording and viewing digital video-EEG over the network, version CLINI C (on request)
- 2.8 Software option for recording and viewing digital video-EEG over a network up to 64 channels, version CLINI C (on request)
- 2.9 Software option for logging and viewing over the network in real time (as required)
- 2.10 Software option for EEG analysis in "on-line" and "off-line" modes, version: full / CFM / DSA. (on demand)
- 2.11 Software option for transmitting EEG data via TCP/IP (optional)
- 2.12 Software option " ERP -manager" (optional)
- 2.13 Software option "Management of outpatient registrar" (as required)
- 3. Software options for viewing:
 - 3.1. EEG software viewing option, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
 - 3.2 Software option for viewing EEG, version: EASY (on request)
 - 3.3 Software option for viewing EEG in "off-line" mode, version: full, CFM, DSA (on request)
 - 3.4 Software option for connecting to databases of external applications (optional)
 - 3.5 EEG software option "Event identification" (as required)
 - 3.6 EEG software option "Analysis of high-frequency oscillations" (on request)
 - 3.7. Software option EEG "Brain Mapping", version: CLINIC, STUDIO (on request)
 - 3.8. EEG software option "Data exchange with MRI" (on request)
 - 3.9. Software EEG option "Interface HOLBERG SCORE" (optional)
 - 3.10. Software EEG option "Audio playback" (optional)
 - 3.11. Software EEG option "Averaging in "off-line" mode" (on demand)
 - 3.12. Software option for viewing digital video-EEG, version: CLINIC, STUDIO (on request)
 - 3.13. Software option for viewing digital video-EEG in Full HD format (on request)

- 3.14. Software option for viewing digital video-EEG in multichannel mode, version: CLINIC (on request)
 - 3.15. Software option for advanced management of the archive server (optional)
 - 3.16. Program statistics option (optional)
 - 3.17. Software option for sleep program interface (optional)
 - 3.18. Software option for background data transfer (optional)
 - 3.19. Software option for exporting EEG data to the Viewer (optional)
 - 3.20. Software option to enable HIS interface (optional)
4. Software options for sleep analysis:
- 4.1 Software option for analyzing imported curves in EDF format (optional)
 - 4.2 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring SD LTM 32 EXPRESS (32 channels) (on request)
 - 4.3 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring SD LTM 64 EXPRESS (64 channels) (optional)
 - 4.4 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring 32 channels (on demand)
 - 4.5 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring 64 channels (on request)
 - 4.6 Software option for the analysis of curves recorded on the EEG Head Units of the SD PLUS series (on request)
5. Program control options.
- 5.1 Software option "Centralized control of video cameras" on a USB -key (on request)
 - 5.2. Software option "Event management" (on demand) on a USB key
6. USB key for registration and viewing programs (optional)
7. USB dongle for viewers for 1/5/10/25 users (optional)
8. USB dongle for sleep analysis programs (optional)
- 9 CB TERMINAL 15 cable for connecting the interface board to the BOX TERMINAL EXPRESS interface module
- EXPRESS wall panel WALL CONNECTOR /12 with additional socket and IN-trigger and OUT-trigger for use in combination with BOX TERMINAL EXPRESS interface module (optional) and 12V voltage amplifier (optional)
- 11 LCD: 19"/22"/24"/27" (optional)
- 12 EEG head unit for long-term monitoring SD LTM 32 EXPRESS (32 channels)
- 13 Electrode box JB LTM 32 EXPRESS to the EEG head unit for long-term monitoring SD LTM 32 EXPRESS (32 channels)
- 14 JB LTM DSUB electrode box with double 37-pin D-Sub connector for EEG headunits for long-term monitoring SD LTM EXPRESS series (on request)
- 15 Electrode box JB LTM TOPHAT (universal) with sensor inputs for all channels (as required)
- 16 LED photostimulator with bracket
- 17 BQ NET POE 2 interface module for connecting SD LTM EXPRESS long-term monitoring EEG headunit. includes SUPP BQ NET POE 2 wall mount and BQ USB EXPRESS interface module

- EXPRESS wall panel WALL CONNECTOR / DE with additional socket and IN-trigger and OUT-trigger for use in combination with BOX TERMINAL EXPRESS interface module (optional)
- 19 BQ NET POWER POE power supply for BQ NET POE series interface modules (optional)
- 20 BQ PCI EXPR LTM FULL interface board based on PCI Express FULL Size to control EEG headunit for long-term monitoring SD LTM EXPRESS (optional)
- 21 BQ PCI EXPR LTM HALF interface board based on PCI Express HALF Size to control the EEG headunit for long-term monitoring SD LTM EXPRESS (optional)
- 22 BQ USB EXPRESS interface module to EEG head unit for long-term monitoring SD LTM EXPRESS (on request)
- 23 BQ USB EXPRESS interface module to EEG head unit for long-term monitoring SD LTM EXPRESS (optional)
- 24 BQ USB PLUS interface module to EEG head unit for long-term monitoring SD LTM EXPRESS (on request)
- 25 USB adapter to serial port ADAPT USB SER (optional)
- 26 Isolated USB to Serial Adapter ADAPT USB SER ISO (as needed)
- 27 CB LTM DC cable for connecting DC channels to EEG Head Units for long-term SD monitoring LTM
- 28 CB MAG RAPID cable for connecting the Micromed system to Magstim magnetic stimulation devices (D-Sub 26 pin (optional)
- 29 Trigger cable CB RCRC4N for stimulator (optional)
- 30 Trigger cable CB RCRC4R for stimulator (optional)
- 31 CB RCRC1R trigger cable for stimulator (optional)
- 32 Pair of CB RCRC5 PAIR trigger cables for stimulator (optional)
- 33 CB SER RCA EXT cable for synchronization with external stimulators (optional)
- 34 Fiber optic cable 10m/15m/20m/30m/40m/50m/100m (optional)
- 35 HBC EXPRESS cable for connecting an interface board or interface module to EEG Head Units (optional)
- 36 Cable HBC INT LTM R for connecting an interface board or interface module to EEG Head Units for long-term monitoring of the SD LTM series (on request)
- 37 Power supplies LED POWER for LED photostimulator
- 38 TERM LTM Quick Disconnect Cable to EEG Head Units
- 39 CAMERA IP HD PRO (optional)
- 40 CAMERA IID Full HD (1920 x 1080) with PTZ control and automatic IR switching (on demand) (on demand)
- 41 CAMERA IP HD OUTDOOR (optional)
- 42 CAMERA IP HD PLUS (optional)
- 43 CAMERA LAPTOP with USB 2.0. including microphone (optional)
- 44 CAMERA IP COMPACT video camera (optional)

- 45 GC MM video card with dual monitor outputs (optional)
- 46 HDMI EXTENDER device for transmitting HDMI signal over Cat 5 e /6 cable, including receiver and transmitter (optional)
- 47 IP MIC for IP cameras (optional)
- 48 IP POWER power adapter for IP cameras (optional)
- 49 VLAMP IR auto-on (on demand)
- 50 Composite switcher VMIXER PRO with 2 inputs (optional)
- 51 Professional switcher VMIXER PRO with 5 HDMI inputs (optional)
- 52 Bracket ARM GM LTM for EEG Head Units for long-term monitoring of the SD LTM series, complete (on request)
- 53 Vertical bracket ARM LCD DOUBLE for two LCD displays for fixing on the cart (optional)
- 54 Horizontal bracket ARM LCD HORIZ for two LCD displays for fixing on the cart (optional)
- 55 ARM LCD SW Vertical LCD Bracket for Trolley Fixing (Optional)
- 56 Adjustable BAG BELT HP for fixing the headunit bag to the patient
- 57 EXPRESS WALL CONNECTOR NT /12 wall panel with additional socket for use in combination with BOX TERMINAL interface module (optional) EXPRESS and 12V voltage booster (optional)
- 58 BAG LTM HP bag for EEG Head Unit for long-term monitoring of the SD LTM series (on request)
- 59 BASE STAND (optional)
- 60 BQ CLAMP (optional)
- 61 CAM JOIN swivel mount for video camera on a stand (optional)
- 62 CAM STAND for mounting a video camera or head unit on a trolley (optional)
- 63 Fixing FIX GM for trolley bracket (optional)
- 64 Fixing FIX GM TW for table bracket (optional)
- 65 JOIN GM connection for attaching the head unit holder to the bracket (optional)
- 66 LTM SCREW screw for attaching the electrode box to the head unit (optional)
- 67 LED STAND PHOTIC wheel base for LED photostimulator (on request)
- 68 LED SUPP PHOTIC bracket for fixing the LED photostimulator (optional)
- 69 STAND CAM SW video camera holder assembly (optional)
- 70 STAND POLE NEW on wheels with weight (optional)
- 71 Wall mount SUPP BQ NET POE 2 wall mount for BQ NET POE 2 interface modules BQ NET POE 2 / BQ NET POE M (optional) (optional) (optional)
- 72 SUPP LTM Holder for EEG Head Units for Long-Term Monitoring of the SD LTM Series (on request)
- 73 VCLAMP mount for video camera or IR light (optional)

- 75 VSUPP wall bracket for video camera (optional)
- 76 VSUPP IP metal support for attaching the IP video camera to the swivel mount (optional)
- 77 VSUPP IP IR metal support for attaching the IP camera to the rack (optional)
- 78 VSUPP IP HD PLUS BLACK wall mount for video camera, black (optional)
- 79 VSUPP IP HD PLUS WHITE wall mount for video camera, white (optional)
- 80 ADAPT BNC RCA adapter (optional)
- 81 ADAPT DVI HDMI adapter (optional)
- 82 ANY WHERE 2 USB over LAN device with 2 USB ports (optional) (optional)
- 83 ANY WHERE 5 USB over LAN device with 5 USB ports (optional)
- 84 USB CAM JOY STICK to control cameras (optional)
- 85 24-key CAM PAD numeric keypad with USB connection for camera control (optional)
- 86 Numeric keypad CAM PAD 48 by 48 keys with USB connection for camera control (optional)
- 87 Connecting cable CB BNBN1 (BNC to BNC) (optional)
- 88 Connecting cable CB BNRC1 (BNC - RCA (Red) (optional)
- 89 Connecting cable CB BNRC5 (BNC - RCA (Black) (optional)
- 90 CB DVI monitor cable (optional)
- 91 Cable HDMI - HDMI MINI for video transmission, length 2m/5m/10m (optional)
- 92 HDMI - HDMI cable for high-definition video transmission, length 2m/5m/15m (optional)
- 93 Audio cable CB JACK JMINI with 6.3 mm plugs - MiniJack 3.5 mm (optional)
- 94 Audio cable CB JACK RCRC with 6.3 mm plugs - double RCA (optional)
- 95 CB MICROPHONE cable (optional)
- 96 CB MON cable for monitor (D-Sub 15 pin - D-Sub 15 pin) (optional)
- 97 CB RCRC JMINI, RCA male to MiniJack 3.5mm two-wire audio cable (optional)
- 98 Two-wire audio cable CB RCRC PAIR, RCA male-double RCA (as required)
- 100 CB VIDEO video cable as required)
- 101 CB IP HD PLUS 15 cable for controlling IP HD PLUS cameras (optional)
- 102 CB IP HD PLUS DC 5 cable for connecting IP HD PLUS cameras to each other (optional)
- 103 CB IP HD PLUS EXTENDER adapter for IP HD PLUS signal transmission over Cat5e/6 cable, pair (optional)
- 104 HDMI BOOSTER self-powered HDMI amplifier (optional)
- 105 CARSTAR mobile trolley (on request)
- 106 FIX TOR bracket for attaching the isolation transformer to the trolley (optional)
- 107 Mouse pad small/regular/large (optional)
- 108 Trolley CAR PLUS, mobile, adjustable in height (as required)

- 109 FIX CART plates for fixing trolley legs (optional)
- 110 Power cable kit KIT CB POWER for trolley (optional)
- 111 KIT LABEL PLUS for trolley (optional)
- 112 SUPP KEYB PLUS keyboard drawer for trolley (optional)
- 113 SUPP PC PLUS trolley mount (optional)
- 114 SUPP PC PLUS SFF small computer trolley mount (optional)
- 115 Horizontal bar BAR POLE NEW for attaching the bracket to the trolley (optional)
- 116 CAR TORO trolley with shelves, standard size (on request)
- 117 CAR TORO SMALL trolley with shelves, small size (on demand)
- 118 COVER TORO rear cover for standard CARTORO trolley (optional)
- 119 Rear cover COVER TORO SMALL for small cart CAR TORO SMALL (optional)
- 120 Isolation transformer TOR 500, Pmax 500 VA, black (on request)
- 121 Isolation transformer TOR 500 W, Pmax 500 VA, white (on request)
- 122 Isolation transformer TR 50, 100 VA (as required)
- 123 NET ISOLAT network isolator for 10/100 Mbit networks, with RJ45 connectors (optional)
- 124 Network cable NET RJ45 with RJ45 connector, length 3m/5m/10m (optional)
- 125 NET SWITCH 4 4-port network switch (optional)
- 126 NET SWITCH 8 8-port network switch (optional)
- 127 NET SWITCH PRO professional network switch (optional)
- 128 Standard Bluetooth USB interface (optional)
- 129 BS USB READER Compact Flash card reader (optional)
- 130 LTM BUTTON EXT extension for LTM BUTTON TP patient button, length 120 cm (optional)
- LTM BUTTON TP Patient Button for use with SD LTM Series EEG Head Units for Long-Term Monitoring
- 132 PLIERS TP forceps for removing electrodes from the electrode box (optional)
- 133 USB extension cable USB EXT CB (optional)
- 134 USB EXTENDER device for transmitting USB signal over Cat5e/6 cable, including receiver and transmitter (optional)
- 135 BQ CONTACT contact closure control module (optional)
- 136 Patient button with 5 m cable and RCA connector (optional)
- 137 BUTTON BOX for controlling patient buttons (optional)
- 138 Cables for BUTTON BOX (optional)
- 139 Illuminated patient button with 2 m cable and LEMO connector (optional)
- 140 Power supply for BUTTON BOX (optional)
- 141 Cable for connection unit for oximetry sensors, Length 3 m (as required)

- 142 Oximetry Sensor Connection Box for use with EEG Head Units for Long-Term Monitoring SD LTM Series (as required)
- 143 Oximetry sensor "EarClip", adult (optional)
- 144 FingerClip oximetry sensor, adult (optional)
- 145 FingerClip oximetry sensor, child (optional)
- 146 Fixing tape SD OXY RA, adult size (on request)
- 147 Fixing tape SD OXY RC, child size (as required)
- 148 Flexible oximetry probe, adult size (optional)
- 149 Flexible oximetry sensor, child size (optional)
- 150 Large reusable soft oximeter probe (digital readout 12.5 - 25.5 mm) (optional)
- 151 Medium size reusable soft oximeter probe (digital readout 10 - 19 mm) (optional)
- 152 Reusable soft oximetry probe small size (digital readout 7.5 - 12.5 mm) (optional)
- 153 COMP DVI RX KVM extender compact, for DVI + USB HID reception (optional)
- 154 KVM Extender COMP DVI RX AUD compact, for receiving DVI + USB HID + analog audio signal (optional)
- 155 KVM Extender COMP DVI TX compact, for DVI + USB HID transmission (optional)
- 156 COMP DVI TX AUD Compact KVM Extender for DVI + USB HID + analog audio transmission (optional)
- 157 KVM DVI module RX for DVI + USB reception HID (on demand)
- 158 KVM DVI module RX AUD dual for DVI + USB reception HID + analog audio signal (optional)
- 159 KVM DVI module TX for DVI + USB transmission HID (on demand)
- 160 KVM DVI module TX AUD dual for DVI + USB transmission HID + analog audio signal (optional)
- 161 KVM DVI module TX AUD DCAT dual for DVI + USB transmission HID + analog audio signal with two-channel transmission (as required)
- 162 KVM DVI module TX DCAT to DVI + USB transmission HID ,dual transmission (optional)
- 163 KIT KVM module including transmitter and receiver for keyboard, mouse and monitor. (optional)
- 164 KVM module KIT KVM DUAL, including transmitter and receiver, for keyboard, mouse and two monitors (optional)
- 165 KVM - KVM Chassis BODY 2 for 2 DVI modules (optional)
- 166 KVM - KVM Chassis BODY 4 for 4 DVI modules (optional)
- 167 KVM -KVM Chassis BODY 4 R for 4 DVI modules , including rack mounting bracket (optional)
- 168 KVM - KVM Chassis BODY 6 for 6 DVI modules (optional)
- 169 KVM -KVM Chassis BODY 6 R for 6 DVI modules , including rack mounting bracket

- 170 IP KVM Module, KVM IP RX, Receiver (optional)
- 171 KVM IP, KVM IP TX, Transmitter (optional)
- 172 KVM matrix KVM SWITCH for 8/16/32/48/64/80 ports (on demand)
- 173 Additional mount SUPP TOPHAT for Electrode box JB LTM TOPHAT (on request)
- 174 EXPRESS WALL CONNECTOR NT /DE wall panel with additional socket for use in combination with BOX TERMINAL EXPRESS interface module (on request)
- 175 HBC BT/RJ45 extension adapter for BOX TERMINAL EXPRESS interface module (module side) (as required)
- 176 HBC HT/RJ45 extension adapter for BOX TERMINAL EXPRESS interface module (PC side) (optional)
- 177 EXTERNAL BATTERY (optional)
- 178 EXTERNAL CHARGER for external battery (optional)
- 179 HBC EXTERNAL BATTERY cable for connecting an external battery (optional)
- 180 180 Extension cable LED EXT for LED photostimulator (on request)
- 181 TERM EXPRESS quick-disconnect cable to the EEG Head Unit for long-term monitoring SD LTM EXPRESS in a configuration of up to 64 channels (on request)
182. Electrode cap (as needed):
- MRINEO 19 X
 - MRINEO 31 X
 - MRINEO 61 X
183. Cable for EL CFM electrodes (as needed):
- CB CFM050
 - CB CFM150
184. Electrodes with cable for LTM (as required):
- CB EL TPCO 19
 - CB EL TPCO 31
185. Electrode (plate) EL CFM for 2-channel EEG (on demand)
186. Disk electrode with cable (as required):
- EL TPCO
 - EL TPCO 6
 - EL TPCO G06
 - EL TPCO G10
 - EL TPCO S 100
 - EL TPCO S 150
 - EL TPCO S 200
 - EL TPCO S 250
 - EL TPCO SG 100
 - EL TPCO SG 150
 - EL TPCO SG 200
 - EL TPCO SG 250
 - EL TPCO SLT 100
 - EL TPCO SLT 150
 - EL TPCO SLT 200
 - EL TPCO SLT 250
 - EL TPCO SLTG 100
 - EL TPCO SLTG 150

-EL _ TPCO SLTG 250

187. Pediatric disc electrode with cable (as needed):

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188. Ear electrode (clip) EL TPEAR with cable (on demand)

189. EL TPLINK cable for combining two ear electrodes EL TPEAR (on request)

190. Connecting cable for LTM / LTM 256 (on request)

191. CB PUD150 cable for EL PUD electrodes (optional)

192. Disposable electrodes EL PUD (on demand)

193. Electrode monopolar needle (on demand):

- EMG MBOT2530
- EMG MBOT2540
- EMG MBOT3745
- EMG MBOT5050
- EMG MBOT7570

194. Snoring sensor (on demand):

- EPM 1200
- EPM 1250

195. EPM 1312 breathing sensor set, includes sensor (EPM 1387) and 4 tapes in different sizes (EPM 1340 + EPM 1341 + EPM 1342 + EPM 1343) (optional)

196. Pediatric breathing sensor EPM 1390 (on request)

197. Pediatric air flow sensor EPM 1400 (on demand)

GM / HM disposable breath sensor kit , includes: EPM 1415 cable and 10 sensors adult 1426 EPMs / medium size 1427 EPMs / infant size 1428 EPMs) (as required)

199. Air flow sensor EPM 1450 (optional)

200. Set of air flow and pressure sensors EPM 1475 (optional)

201. Air pressure and snoring sensor EPM 14830 (optional)

202. Body position sensor EPM 1563 (optional)

203. EPM 1715 motion sensor kit (including piezo transducer (EPM 1770) and pairs of ankle straps (EPMS 1730) (optional)

204. Chest force sensor kit, including interface cable and 4 chest straps (as required):

- EPM 9102 A for adults
- EPM 9102 C pediatric

205. Abdominal force sensor kit, including interface cable and 4 abdominal bands (as required):

- EPM 9101 A for adults
- EPM 9101 C pediatric

206. A set of abdominal and thoracic force sensors, including abdominal and thoracic interface cables, and 4 tapes (as required):

- EPM 914x A for adults
- EPM 914x C pediatric
- EPM 915x A for adults
- EPM 915x C pediatric

207. Tape for the breath sensor (as needed):

- EPMs 1340
- EPMs 1341
- EPMs 1342
- EPMs 1343

208. Piezocrystalline breath sensor EPMs 1387 (optional)
209. Cable for single-channel airflow sensor EPMs 1415 (optional)
210. Thermal breath sensor (as required):
- EPMs 1426 adult size
 - EPMs 1427 medium size
 - EPMs 1428 pediatric
211. EPMS 1450-B interface for EPM 1450 airflow sensor (optional)
- S transducer for EPM 1450 airflow sensor (optional)
213. Combined thermistor pressure and air flow sensor EPMs 1475
- CO 2 cannula (as required):
- EPMs 14802, nasal and oral, male adult
 - EPMs 14803, male for newborns
 - EPMs 14804, male pediatric
 - EPMs 14805, male adult
215. EPMs 14830 piezoelectric airflow and snore pressure transducer (optional)
216. Elastic tape EPMs 1550 for body position sensor EPM 1520 (optional)
217. EPMs 1563 transducer for EPM 1563 body position sensor (optional)
218. ThermoCan cannula (optional):
- EPMs 15802
 - EPMs 15807
219. Pair of EPMs 1730 ankle straps (optional)
220. Piezoelectric transducer EPMs 1770 for motion sensor EPM 1715 (optional)
221. A pair of induction tapes (as needed):
- 9002-L120 EPMs
 - 9002-L133 EPMs
 - 9002-L40 EPMs
 - 9002-L60 EPMs
 - 9002-L90 EPMs
222. Adjustable induction band EPMs 9010 (optional)
223. Induction interface cable (as required):
- EPMs 9101 for the abdomen
 - EPMs 9102 for the thoracic region
224. Disposable ECG electrode ECG EL (on request)
225. Cable for ECG electrode (optional):
- ECG EL TPCLB black
 - ECG EL TPCLR red

II. ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL STUDIES MICROMED , PERFORMED BY: SD LTM 32 PLUS , CONSISTING OF:

System software PLUS EVOLUTION - BRAIN QUICK including the following software options:

Software option for registration and viewing:

- 2.1 Software option for recording and viewing EEG, 2-channel CFM version (on request)
- 2.2 Software option for recording and viewing EEG, 8-channel CFM version (on request)
- 2.3 Software option for recording and viewing EEG, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
- 2.4 Software option for recording and viewing digital video-EEG up to 64 channels, version CLINI C (on request)
- 2.5 Software option for recording and viewing digital video-EEG, version CLINI C (on request)

- 2.5. Software option for recording and viewing digital video-EEG in Full HD format (on request)
- 2.7 Software option for recording and viewing digital video-EEG over the network, version CLINI C (on request)
- 2.8 Software option for recording and viewing digital video-EEG over a network up to 64 channels, version CLINI C (on request)
- 2.9 Software option for logging and viewing over the network in real time (as required)
- 2.10 Software option for EEG analysis in "on-line" and "off-line" modes, version: full / CFM / DSA. (on demand)
- 2.11 Software option for transmitting EEG data via TCP/IP (optional)
- 2.12 Software option "ERP -manager" (optional)
- 2.13 Software option "Management of outpatient registrar" (as required)
- 3. Software options for viewing:
 - 3.1. EEG software viewing option, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
 - 3.2 Software option for viewing EEG, version: EASY (on request)
 - 3.3 Software option for viewing EEG in "off-line" mode, version: full, CFM, DSA (on request)
 - 3.4 Software option for connecting to databases of external applications (optional)
 - 3.5 EEG software option "Event identification" (as required)
 - 3.6 EEG software option "Analysis of high-frequency oscillations" (on request)
 - 3.7. Software option EEG "Brain Mapping", version: CLINIC, STUDIO (on request)
 - 3.8. EEG software option "Data exchange with MRI" (on request)
 - 3.9. Software EEG option "Interface HOLBERG SCORE" (optional)
 - 3.10. Software EEG option "Audio playback" (optional)
 - 3.11. Software EEG option "Averaging in "off-line" mode" (on demand)
 - 3.12. Software option for viewing digital video-EEG, version: CLINIC, STUDIO (on request)
 - 3.13. Software option for viewing digital video-EEG in Full HD format (on request)
 - 3.14. Software option for viewing digital video-EEG in multichannel mode, version: CLINIC (on request)
 - 3.15. Software option for advanced management of the archive server (optional)
 - 3.16. Program statistics option (optional)
 - 3.17. Software option for sleep program interface (optional)
 - 3.18. Software option for background data transfer (optional)
 - 3.19. Software option for exporting EEG data to the Viewer (optional)
 - 3.20. Software option to enable IIS interface (optional)
- 4. Software options for sleep analysis:
 - 4.1 Software option for analyzing imported curves in EDF format (optional)
 - 4.2 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring 32 channels (on request)
 - 4.3 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring 64 channels (on request)
 - 4.4 Software option for the analysis of curves recorded on the EEG Head Units of the SD PLUS series (on request)
- 5. Program control options.
 - 5.1 Software option "Centralized control of video cameras" on a USB -key (on request)
 - 5.2. Software option "Event management" USB -key (as required)

6. USB key for registration and viewing programs (optional)
7. USB dongle for viewers for 1/5/10/25 users (optional)
8. USB dongle for sleep analysis programs (optional)
- 9 HBC EXPRESS cable for connecting an interface board or interface module to EEG Head Units (optional)
- 10 BQ interface module USB PLUS to EEG headunit for long-term SD monitoring LTM PLUS (64 channels)
- 11 LCD: 19"/22"/24"/27" (optional)
- 12 EEG head unit for long-term SD monitoring LTM 32 PLUS (32 channels)
- 13 Electrode box JB LTM 32 PLUS to the EEG head unit for long-term SD monitoring LTM 32 PLUS (32 channels)
14. Electrode box JBD SUB PLUS with Dual 37-Pin D - Sub Connector to EEG Head Units for SD Series Long-Term Monitoring LTM PLUS
- 15 JB TOPHAT PLUS electrode box for all channels with double 37-pin D-Sub connector for EEG head units for long-term monitoring of the SDLTMPLUS series
- 16 BOX TERMINAL MULTI interface module for connecting the EEG head unit for long-term SD monitoring LTM PLUS (up to 256 channels)
- 17 LED photostimulator with bracket
- 18 BQ NET POE PLUS interface module for connecting EEG head unit for long-term monitoring SD LTM PLUS , includes SUPP wall mount BQ NET POE 2
- 19 BQ NET POWER POE power supply for BQ NET POE series interface modules (optional)
- 20 USB adapter to serial port ADAPT USB SER (optional)
- 21 Isolated USB to Serial Adapter ADAPT USB SER (optional) ISO
- 22 CB LTM DC cable for connecting DC channels to EEG Head Units for long-term monitoring SDLTM
- 23 CB MAG RAPID cable for connecting the Micromed system to Magstim magnetic stimulation devices (D-Sub 26 pin (optional)
- 24 Trigger cable CB RCRC4N for stimulator (optional)
- 25 Trigger cable CB RCRC4R for stimulator (optional)
- 26 CB RCRC1R trigger cable for stimulator (optional)
- 27 Pair of CB RCRC5 PAIR trigger cables for stimulator (optional)
- 28 CB SER RCA EXT cable for synchronization with external stimulators (optional)
- 29 Fiber optic cable 10m/15m/20m/30m/40m/50m/100m (optional)
- 30 Cable HBC INT LTM R to connect the interface board or interface module to EEG Head Units for long-term monitoring of the SDLTM series
- 31 Power supplies LED POWER for LED photostimulator
- 32180 Extension cable LED EXT for LED photostimulator (on request)

- 34 CAMERA IP HD PRO (optional)
- 35 CAMERA HD Full HD (1920 x 1080) with PTZ control and automatic IR switching (on demand) (on demand)
- 36 CAMERA IP HD OUTDOOR (optional)
- 37 CAMERA IP HD PLUS (optional)
- 38 CAMERA LAPTOP with USB 2.0, including microphone (optional)
- 39 CAMERA IP COMPACT video camera (optional)
- 40 GC MM video card with dual monitor outputs (optional)
- 41 HDMI EXTENDER device for transmitting HDMI signal over Cat 5 e /6 cable, including receiver and transmitter (optional)
- 42 IP MIC for IP cameras (optional)
- 43 IP POWER power adapter for IP cameras (optional)
- 44 VLAMP IR auto-on (on demand)
- 45 Composite switcher VMIXER PRO with 2 inputs (optional)
- 46 Professional switcher VMIXER PRO with 5 HDMI inputs (optional)
- 47 Bracket ARM GM LTM for EEG head units for long-term monitoring of the SD LTM series, complete (on request)
- 48 Vertical bracket ARM LCD DOUBLE for two LCD displays for fixing on the cart (optional)
- 49 Horizontal bracket ARM LCD HORIZ for two LCD displays for fixing on the cart (optional)
- 50 ARM LCD SW Vertical LCD Bracket for Trolley Fixing (optional)
- 51 Adjustable BAG BELT HP for fixing the headunit bag to the patient
- 52 TERM PLUS quick-disconnect cable to the EEG headunit for long-term monitoring SD LTM PLUS in a configuration of up to 64 channels
- 53 BAG LTM HP Bag for EEG Head Unit for Long-Term Monitoring SD Series LTM
- 54 BASE STAND (optional)
- 55 BQ CLAMP (optional)
- 56 Rotary mount CAM JOIN for a video camera on a stand (optional)
- 57 CAM STAND for mounting a video camera or head unit on a trolley (optional)
- 58 Fixing FIX GM for trolley bracket (optional)
- 59 Fixing FIX GM TW for table bracket (optional)
- 60 JOIN GM connection for attaching the head unit holder to the bracket (optional)
- 61 LTM SCREW screw for attaching the electrode box to the head unit (optional)
- 62 LED STAND PHOTIC wheel base for LED photostimulator (on request)
- 63 LED SUPP PHOTIC bracket for fixing the LED photostimulator (optional)
- 64 STAND CAM SW video camera holder assembly (optional)

- 66 Wall mount SUPP BQ NET P OE2 for interface modules BQNETPOE 2 / BQ NET POE M (as required)
- 67 SUPP LTM Holder for SD LTM Series EEG Head Units (optional)
- 68 ADAPT DVI HDMI adapter (optional)
- 69 VCLAMP mount for video camera or IR light (optional)
- 70 VSTAND tripod stand for video camera or IR light (optional)
- 71 VSUPP wall bracket for video camera (optional)
- 72 VSUPP IP metal support for attaching the IP video camera to the swivel mount (optional)
- 73 VSUPP IP IR metal support for attaching the IP camera to the rack (optional)
- 74 VSUPP IP HD PLUS BLACK wall mount for video camera, black (optional)
- 75 VSUPP IP HD PLUS WHITE wall mount for video camera, white (optional)
- 76 ADAPT BNC RCA adapter (optional)
- 78 ANY WHERE 2 USB over LAN device with 2 USB ports (optional)
- 79 ANY WHERE 5 USB over LAN device with 5 USB ports (optional)
- 80 USB CAM JOYSTICK to control cameras (optional)
- 81 24-key CAM PAD numeric keypad with USB connection for camera control (optional)
- 82 Numeric keypad CAM PAD 48 by 48 keys with USB connection for camera control (optional)
- 83 Connecting cable CB BNBNI (BNC to BNC) (optional)
- 84 Connecting cable CB BNRC1 (BNC - RCA (Red) (optional)
- 85 Connecting cable CB BNRC5 (BNC - RCA (Black) (optional)
- 86 CB DVI monitor cable (optional)
- 87 Cable HDMI - HDMI MINI for video transmission, length 2m/5m/10m (optional)
- 88 HDMI to HDMI cable for high-definition video transmission, 2m/5m/15m (optional)
- 89 Audio cable CB JACK JMINI with 6.3 mm plugs - MiniJack 3.5 mm (optional)
- 90 Audio cable CB JACK RCRC with 6.3 mm plugs - double RCA (optional)
- 91 CB MICROPHONE cable (optional)
- 92 CB MON cable for monitor (D-Sub 15 pin - D-Sub 15 pin) (optional)
- 93 CB RCRC JMINI, RCA male to MiniJack 3.5mm two-wire audio cable (optional)
- 94 Two-wire audio cable CB RCRC PAIR, RCA male to double RCA (as required)
- 95 CB VIDEO video cable (optional)
- 96 CB IP HD PLUS 15 cable for controlling IP cameras (optional)
- 97 CB IP HD PLUS DC 5 cable for connecting IP HD PLUS cameras to each other (optional)
- 98 CB IP HD PLUS EXTENDER adapter for IP HD PLUS signal transmission over Cat5e/6 cable, pair (optional)

- 100 CARSTAR mobile trolley (on request)
- 101 FIX TOR bracket for attaching the isolation transformer to the trolley (optional)
- 102 Mouse pad small/regular/large (optional)
- 103 CAR PLUS trolley, mobile, height-adjustable (on demand)
- 104 FIX CART plates for fixing trolley legs (optional)
- 105 Power cable kit KIT CB POWER for trolley (optional)
- 106 KIT LABEL PLUS for trolley (optional)
- 107 SUPP KEYB PLUS keyboard drawer for trolley (optional)
- 108 SUPP PC PLUS trolley mount (optional)
- 109 Small computer mount SUPP PC PLUS SFF on cart (optional)
- 110 Horizontal bar BAR POLE NEW for attaching the bracket to the trolley (optional)
- 111 CAR TORO trolley with shelves, standard size (on request)
- 112 CAR TORO SMALL trolley with shelves, small size (on request)
- 113 COVER TORO rear cover for standard CARTORO trolley (optional)
- 114 Rear cover COVER TORO SMALL for small cart CAR TORO SMALL (optional)
- 115 Isolation transformer TOR 500, Pmax 500 VA, black (on request)
- 116 Isolation transformer TOR 500 W, Pmax 500 VA, white (on request)
- 117 Isolation transformer TR 50, 100 VA (as required)
- 118 NET ISOLAT network isolator for 10/100 Mbit networks, with RJ45 connectors (optional)
- 119 Network cable NET RJ45 with RJ45 connector, length 3m/5m/10m (optional)
- 120 NET SWITCH 4 4-port network switch (optional)
- 121 NET SWITCH 8 8-port network switch (optional)
- 122 NET SWITCH PRO professional network switch (optional)
- 123 Standard Bluetooth USB interface (optional)
- 124 BS USB READER Compact Flash card reader (optional)
- 125 Extension LTM BUTTON EXT for patient button LTM BUTTON TP, length 120 cm (optional)
- LTM BUTTON TP Patient Button for use with SD LTM Series EEG Head Units for Long-Term Monitoring
- 127 PLIERS TP forceps for removing electrodes from the electrode box (optional)
- 128 USB extension cable USB EXT CB (optional)
- 129 USB EXTENDER device for transmitting USB signal over Cat5e/6 cable, including receiver and transmitter (optional)
- 130 BQ CONTACT contact closure control module (optional)
- 131 Patient button with 5 m cable and RCA connector (optional)

- 133 Cables for BUTTON BOX (optional)
- 134 Illuminated patient button with 2 m cable and LEMO connector (optional)
- 135 Power supply for BUTTON BOX (optional)
- 136 Cable for oximetry sensor connection unit. Length 3 m (as required)
- 137 Oximetry Sensor Connection Box for use with EEG Head Units for Long-Term Monitoring SD Series LTM (on demand)
- 138 Oximetry sensor "EarClip", adult (optional)
- 139 FingerClip oximetry sensor, adult (optional)
- 140 FingerClip oximetry sensor, child (optional)
- 141 Fixing tape for SD OXY RA, adult size (on request)
- 142 Fixing tape SD OXY RC, child size (as required)
- 143 Flexible oximetry probe, adult size (optional)
- 144 Flexible oximetry sensor, child size (optional)
- 145 Large reusable soft oximeter probe (digital readout 12.5 - 25.5 mm) (optional)
- 146 Medium size reusable soft oximeter probe (digital readout 10 - 19 mm) (optional)
- 147 Reusable soft oximeter probe small size (digital readout 7.5 - 12.5 mm) (optional)
- 148 COMP DVI RX KVM extender compact, for DVI + USB HDI reception (optional)
- 149 KVM Extender COMP DVI RX AUD compact, for receiving DVI + USB HDI + analog audio signal (optional)
- 150 KVM Extender COMP DVI TX compact, for DVI + USB HDI transmission (optional)
- 151 COMP DVI TX AUD Compact KVM Extender for DVI + USB HDI + analog audio transmission (optional)
- 152 KVM DVI module RX for DVI + USB reception HID (on demand)
- 153 KVM DVI module RX AUD dual for DVI + USB reception HID + analog audio signal (optional)
- 154 KVM DVI module TX for DVI + USB transmission HID (on demand)
- 155 KVM DVI module TX AUD dual for DVI + USB transmission HID + analog audio signal (optional)
- 156 KVM DVI module TX AUD DCAT dual for DVI + USB transmission HID + analog audio signal with two-channel transmission (as required)
- 157 KVM - DVI TX DCAT Module for DVI + USB HID Dual Link Transmission
- 158 KIT KVM module including transmitter and receiver for keyboard, mouse and monitor. (optional)
- 159 KVM module KIT KVM DUAL, including transmitter and receiver, for keyboard, mouse and two monitors (optional)
- 160 KVM -KVM Chassis BODY 2 for 2 DVI modules (optional)
- 161 KVM -KVM Chassis BODY 4 for 4 DVI modules (optional)

162 KVM -KVM Chassis BODY 4 R for 4 DVI modules , including rack mounting bracket (optional)

163 KVM - KVM Chassis BODY 6 for 6 DVI modules (optional)

164 KVM - KVM Chassis BODY 6 R for 6 DVI modules , including rack mounting bracket (optional)

165 IP KVM Module. KVM IP RX, Receiver (optional)

166 KVM IP, KVM IP TX, Transmitter (optional)

167 KVM matrix KVM SWITCH for 8/16/32/48/64/80 ports (on demand)

168.External battery EXTERNAL BATTERY(on demand)

169. EXTERNAL CHARGER charger for external battery (on request)

170.HBC EXTERNAL BATTERY cable for connecting an external battery (optional)

171. Electrode cap (as needed):

- MRINEO 19 X

- MRINEO 31 X

- MRINEO 61 X

172. Cable for EL CFM electrodes (as required):

- CB CFM050

- CB CFM150

173. Electrodes with cable for LTM (as required):

- CB EL TPCO 19

- CB EL TPCO 31

174. Electrode (plate) EL CFM for 2-channel EEG (on demand)

175. Disk electrode with cable (as needed):

- EL TPCO

- EL TPCO 6

- EL TPCO G06

- EL TPCO G10

- EL TPCO S 100

- EL TPCO S 150

- EL TPCO S 200

- EL TPCO S 250

- EL TPCO SG 100

- EL TPCO SG 150

- EL TPCO SG 200

- EL TPCO SG 250

- EL TPCO SLT 100

- EL TPCO SLT 150

- EL TPCO SLT 200

- EL TPCO SLT 250

- EL TPCO SLTG 100

- EL TPCO SLTG 150

- EL TPCO SLTG 200

- EL _ TPCO SLTG 250

176. Pediatric disk electrode with cable (as required):

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- 177. Ear electrode (clip) EL TPEAR with cable (on request)
- 178. EL TPLINK cable for combining two ear electrodes EL TPEAR (optional)
- 179. Connecting cable for LTM / LTM 256 (on request)
- 180. CB PUD150 cable for EL PUD electrodes (optional)
- 181. Disposable electrodes EL PUD (on demand)
- 182. Electrode monopolar needle (on demand):

- EMG MBOT2530
- EMG MBOT2540
- EMG MBOT3745
- EMG MBOT5050
- EMG MBOT7570

- 183. Snoring sensor (on demand):

- EPM 1200
- EPM 1250

- 184. EPM 1312 breathing sensor set, includes sensor (EPM 1387) and 4 tapes in different sizes (EPM 1340 + EPM 1341 + EPM 1342 + EPM 1343) (optional)

- 185. Pediatric breathing sensor EPM 1390 (on request)

- 186. Pediatric air flow sensor EPM 1400 (on demand)

- 187. EPM 1439FM/ GM / HM Disposable Breath Sensor Set, includes: EPM 1415 cable and 10 sensors for Adult 1426 EPMs / Medium Size 1427 EPMs / Infant Size 1428 EPMs) (as required)

- 188. Air flow sensor EPM 1450 (optional)

- 189. Set of air flow and pressure sensors EPM 1475 (optional)

- 190. Air pressure and snoring sensor EPM 14830 (optional)

- 191. Body position sensor EPM 1563 (optional)

- 192. EPM 1715 motion sensor kit (including piezo transducer (EPM 1770) and pairs of ankle straps (EPMS 1730) (optional)

- 193. Chest force sensor kit, including interface cable and 4 chest straps (as required):

- EPM 9102 A for adults
- EPM 9102 C pediatric

- 194. Abdominal force sensor kit, including interface cable and 4 abdominal bands (as required):

- EPM 9101 A for adults
- EPM 9101 C pediatric

- 195. A set of abdominal and thoracic force sensors, including abdominal and thoracic interface cables, and 4 tapes (as required):

- EPM 914x A for adults
- EPM 914x C pediatric
- EPM 915x A for adults
- EPM 915x C pediatric

- 196. Tape for the breath sensor (as needed):

- EPMs 1340
- EPMs 1341
- EPMs 1342
- EPMs 1343

- 197. Piezocrystalline breath sensor EPMs 1387 (optional)

- 198. Cable for single-channel airflow sensor EPMs 1415 (optional)

- 199. Thermal breath sensor (as needed):

- EPMs 1426 adult size
- EPMs 1427 medium size
- EPMs 1428 pediatric

- 200. EPMs 1450-B interface for EPM 1450 airflow sensor (optional)

- 23
- S transducer for EPM 1450 air flow sensor (optional)
202. Combined thermistor pressure and air flow sensor EPMs 1475
CO 2 cannula (as required):
- EPMs 14802, nasal and oral, male adult
 - EPMs 14803, male for newborns
 - EPMs 14804, male pediatric
 - EPMs 14805, male adult
204. EPMs 14830 piezoelectric airflow and snore pressure transducer (optional)
205. Elastic tape EPMs 1550 for body position sensor EPM 1520 (optional)
206. EPMs 1563 transducer for EPM 1563 body position sensor (optional)
207. ThermoCan cannula (optional):
- EPMs 15802
 - EPMs 15807
208. Pair of EPMs 1730 ankle straps (optional)
209. Piezoelectric transducer EPMs 1770 for motion sensor EPM 1715 (optional)
210. A pair of induction tapes (as needed):
- 9002-L120 EPMs
 - 9002-L133 EPMs
 - 9002-L40 EPMs
 - 9002-L60 EPMs
 - 9002-L90 EPMs
211. Adjustable induction tape EPMs 9010 (as required)
212. Induction interface cable (as required):
- EPMs 9101 for the abdomen
 - EPMs 9102 for the thoracic region
213. Disposable ECG electrode ECG EL (on request)
214. Cable for ECG electrode (optional):
- ECG EL TPCLB black
 - ECG EL TPCLR red

III. ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL STUDIES MICROMED , IN VERSIONS: SD LTM 64 EXPRESS , CONSISTING OF:

System software PLUS EVOLUTION - BRAIN QUICK including the following software options:

Software option for registration and viewing:

- 2.1 Software option for recording and viewing EEG, 2-channel CFM version (on request)
- 2.2 Software option for recording and viewing EEG, 8-channel CFM version (on request)
- 2.3 Software option for recording and viewing EEG, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
- 2.4. Software option for recording and viewing digital video-EEG up to 128 channels, version: CLINI C (on request)
- 2.5. Software option for recording and viewing digital video-EEG up to 256 channels, CLINI C version (on request)
- 2.6 Software option for recording and viewing digital video-EEG up to 64 channels, version CLINI C (on request)
- 2.7 Software option for recording and viewing digital video-EEG, version CLINI C (on request)
- 2.8. Software option for recording and viewing digital video-EEG in Full HD format (on request)

request)

- 2.9 Software option for recording and viewing digital video-EEG over the network, version CLINIC (on request)
- 2.10 Software option for recording and viewing digital video-EEG over a network of up to 128 channels, version CLINIC (on request)
- 2.11 Software option for recording and viewing digital video-EEG over a network up to 256 channels, CLINIC version (on request)
- 2.12 Software option for recording and viewing digital video-EEG over a network of up to 64 channels, version CLINIC (on request)
- 2.13 Software option for real-time logging and viewing over the network (as required)
- 2.14 Software option for EEG analysis in "on-line" and "off-line" modes. version: full / CFM / DSA. (on demand)
- 2.15 Software option for transmitting EEG data via TCP/IP (optional)
- 2.16 Software option "ERP -manager" (optional)
- 2.17 Software option "Management of outpatient registrar" (as required)
- 3. Software options for viewing:
 - 3.1. EEG software viewing option, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
 - 3.2 Software option for viewing EEG, version: EASY (on request)
 - 3.3 Software option for viewing EEG in "off-line" mode, version: full, CFM. DSA (on request)
 - 3.4 Software option for connecting to databases of external applications (optional)
 - 3.5 EEG software option "Event identification" (as required)
 - 3.6 EEG software option "Analysis of high-frequency oscillations" (on request)
 - 3.7. Software option EEG "Brain Mapping", version: CLINIC, STUDIO (on request)
 - 3.8. EEG software option "Data exchange with MRI" (on request)
 - 3.9. Software EEG option "Interface HOLBERG SCORE" (optional)
 - 3.10. Software EEG option "Audio playback" (optional)
 - 3.11. Software EEG option "Averaging in "off-line" mode" (on demand)
 - 3.12. Software option for viewing digital video-EEG, version: CLINIC, STUDIO (on request)
 - 3.13. Software option for viewing digital video-EEG in Full HD format (on request)
 - 3.14. Software option for viewing digital video-EEG in multichannel mode, version: CLINIC (on request)
 - 3.15. Software option for advanced management of the archive server (optional)
 - 3.16. Program statistics option (optional)
 - 3.17. Software option for sleep program interface (optional)
 - 3.18. Software option for background data transfer (optional)
 - 3.19. Software option for exporting EEG data to the Viewer (optional)
 - 3.20. Software option to enable IIS interface (optional)
- 4. Software options for sleep analysis:
 - 4.1 Software option for analyzing imported curves in EDF format (optional)
 - 4.2 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring SD LTM 32 EXPRESS (32 channels) (on request)
 - 4.3 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring SD LTM 64 EXPRESS (64 channels) (optional)
 - 4.4 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring 32 channels (on demand)
 - 4.5 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring 64 channels (on request)

5. Program control options.

5.1 Software option "Centralized control of video cameras" on a USB -key (on request)

5.2 Software option "Event management" on a USB key (optional)

6. USB key for registration and viewing programs (optional)

7. USB dongle for viewers for 1/5/10/25 users (optional)

8. USB dongle for sleep analysis programs (optional)

9 EXTERNAL BATTERY (optional)

10 EXTERNAL CHARGER for external battery (optional)

11 LCD: 19"/22"/24"/27" (optional)

12 EEG head unit for long-term monitoring SD LTM 64 EXPRESS (64 channels)

13 Electrode box JB LTM 64 EXPRESS to the EEG head unit for long-term monitoring SD LTM 64 EXPRESS (64 channels)

14 JB LTM DSUB electrode box with double 37-pin D-Sub connector for EEG headunits for long-term monitoring SD LTM EXPRESS series (on request)

15 Electrode box JB LTM TOPHAT (universal) with sensor inputs for all channels (as required)

16 SD LTM STIM Cortical Stimulator to EEG Head Units for Long-Term Monitoring SD Series LTM (64 channels)

17 LED photostimulator with bracket

18 BQ NET POE 2 interface module for connecting SD LTM EXPRESS long-term monitoring EEG headunit. includes SUPP BQ NET POE 2 wall mount and module

BQ power supply NET POWER POE for BQ series interface modules NET POE

20 BQ PCI EXPR LTM FULL interface board based on PCI Express FULL Size to control EEG headunit for long-term monitoring SD LTM EXPRESS (optional)

BQ interface board PCI EXPR LTM HALF , based on PCI Express HALF Size to control EEG head unit for long-term SD monitoring LTM EXPRESS

BQ interface module USB EXPRESS to EEG Headunit SD PLUS and EEG Headunit for long-term SD monitoring LTM EXPRESS

23 BQ USB EXPRESS interface module to EEG head unit for long-term monitoring SD LTM EXPRESS (optional)

BQ interface module USB PLUS to EEG headunit for long-term SD monitoring LTM EXPRESS

25 Voltage adapter ADAPT LTM BOX to EEG head units for long-term monitoring SD LTM EXPRESS (128 and 256-channel systems) connected to the BOX interface module TERMINAL EXPRESS

26 External voltage adapter ADAPT LTM HTHT to SD LTM EXPRESS long-term monitoring EEG headunits (128 and 256-channel systems) connected directly to the BQ PCI EXPR LTM interface board (optional)

27 USB adapter to serial port ADAPT USB SER (optional)

28 Isolated USB to Serial Adapter ADAPT USB SER (optional) ISO

- 29 CB LTM cable for connecting SD cortical stimulator LTM STIM to EEG head unit for long-term SD monitoring LTM64 _
- 30 CB LTM DC cable for connecting DC channels to EEG Head Units for long-term SD monitoring LTM (on demand)
- 31 CB LTM SYNC cable for connecting two EEG Head Units for long-term SD monitoring LTM 64 EXPRESS (optional)
- 32 CB MAG RAPID cable for connecting the Micromed system to Magstim magnetic stimulation devices (D-Sub 26 pin (optional)
- 33 CB RCRC4N trigger cable for stimulator (optional)
- 34 Trigger cable CB RCRC4R for stimulator (optional)
- 35 CB RCRC1R trigger cable for stimulator (optional)
- 36 Pair of CB RCRC5 PAIR trigger cables for stimulator (optional)
- 37 CB SER RCA EXT cable for synchronization with external stimulators
- 38 CB TERMINAL 15 cable for connecting the interface board to the BOX TERMINAL EXPRESS interface module
- 39 EXPRESS wall panel WALL CONNECTOR /12 with additional socket and IN-trigger and OUT-trigger for use in combination with BOX TERMINAL EXPRESS interface module (optional) and 12V voltage amplifier (optional)
- EXPRESS wall panel WALL CONNECTOR / DE with additional socket and IN-trigger and OUT-trigger for use in combination with BOX TERMINAL EXPRESS interface module (optional)
- 41 EXPRESS WALL CONNECTOR NT /12 wall panel with additional socket for use in combination with BOX TERMINAL interface module (optional)EXPRESS and 12V voltage booster (optional)
- 42EXPRESS WALL CONNECTOR NT /DE wall panel with additional socket for use in combination with the BOX TERMINAL EXPRESS interface module (on request)
- 43 Fiber optic cable 10m/15m/20m/30m/40m/50m/100m (optional)
- 44 HBC BT/RJ45 extension adapter for BOX TERMINAL EXPRESS interface module (module side) (optional)
- 45 HBC HT/RJ45 extension adapter for BOX TERMINAL EXPRESS interface module (PC side) (optional)
- 46 HBC EXPRESS cable for connecting an interface board or interface module to EEG Head Units (optional)
- 47 HBC INT LTM R cable for connecting an interface board or interface module to EEG Head Units for long-term monitoring SD series LTM
- 48 Voltage adapter LTU (12V) for mounting inside the BOX TERMINAL EXPRESS interface module to the EEG head unit for long-term monitoring SD LTM EXPRESS in a configuration for 128 and 256 channels
- 49 External voltage adapter LTU HT (12 V) to the EEG Head Unit for long-term monitoring SD LTM in 128 and 256 channel configurations
- 50 LED POWER power supplies for LED photostimulator

- 51 Extension cable LED EXT for LED photostimulator (on request)
- 52 TERM EXPRESS quick-disconnect cable to the EEG headunit for long-term monitoring SD LTM EXPRESS in a configuration of up to 64 channels (on request)
- 53 TERM EXPRESS 128 quick-disconnect cable to two EEG Head Units for long-term monitoring SD LTM EXPRESS in a 128-channel configuration
- 54 Quick disconnect cable TERM EXPRESS 256 to four EEG Head Units for long-term monitoring SD LTM EXPRESS in a 256-channel configuration
- 55 TERM LTM Quick Disconnect Cable to EEG Head Units
- 56 CAMERA IP HD PRO (optional)
- 57 CAMERA HD Full HD (1920 x 1080) with PTZ control and automatic IR switching (on demand) (on demand)
- 58 CAMERA IP HD OUTDOOR (optional)
- 59 CAMERA IP HD PLUS (optional)
- 60 CAMERA LAPTOP with USB 2.0, including microphone (optional)
- 61 CAMERA IP COMPACT video camera (optional)
- 62 GC MM video card with dual monitor output (optional)
- 63 HDMI EXTENDER device for transmitting HDMI signal over Cat 5 e /6 cable, including receiver and transmitter (optional)
- 64 IP MIC for IP cameras (optional)
- 65 IP POWER power adapter for IP cameras (optional)
- 66 VLAMP IR auto-on (on demand)
- 67 VMIXER PRO 2-input composite switcher (optional)
- 68 VMIXER PRO professional switcher with 5 HDMI inputs (optional)
- 69 Bracket ARM GM LTM for EEG head units for long-term monitoring of the SD LTM series, assembled (on request)
- 70 Vertical bracket ARM LCD DOUBLE for two LCD displays for fixing on the cart (optional)
- 71 Horizontal bracket ARM LCD HORIZ for two LCD displays for fixing on the cart (optional)
- 72 ARM LCD SW Vertical LCD Bracket for Trolley Fixing (Optional)
- 73 Adjustable BAG BELT HP for fixing the headunit bag to the patient
- 74 Bag BAG LTM 128 for two EEG head units for long-term monitoring of the SD LTM series (on request)
- 75 BAG LTM HP Bag for EEG Head Unit for Long-Term Monitoring SD Series LTM
- 76 BASE STAND (optional)
- 77 BQ CLAMP (optional)
- 78 CAM JOIN swivel mount for video camera on a stand (optional)
- 79 CAM STAND for mounting a video camera or head unit on a trolley (optional)

- 81 Fixture FIX GM TW for table bracket (optional)
- 82 JOIN GM connection for attaching the head unit holder to the bracket (optional)
- 83 LTM SCREW screw for attaching the electrode box to the head unit (optional)
- 84 LED STAND PHOTIC wheelbase for LED photostimulator (optional)
- 85 LED SUPP PHOTIC bracket for fixing the LED photostimulator (optional)
- 86 STAND CAM SW video camera holder assembly (optional)
- 87 STAND POLE NEW on wheels with weight (optional)
- 88 SUPP BQ NET POE2 wall mount for BQ interface modules NET POE2 / BQ NET POE M (as needed)
- 89 SUPP LTM Holder for EEG Head Units for Long-Term Monitoring of the SD LTM Series (optional)
- 90 Holder SUPP LTM 128 for two EEG Head Units for long-term monitoring of the SD LTM series in a 128-channel configuration
- 91 Holder SUPP LTM 256 for four EEG Head Units for long-term monitoring of the SD LTM series in a 256-channel configuration
- 92 VCLAMP mount for video camera or IR light (optional)
- 93 VSTAND tripod stand for video camera or IR light (optional)
- 94 VSUPP wall bracket for video camera (optional)
- 95 VSUPP IP metal support for attaching the IP video camera to the swivel mount (optional)
- 96 VSUPP IP IR metal support for attaching the IP camera to the rack (optional)
- 97 VSUPP IP HD PLUS BLACK wall mount for video camera, black (optional)
- 98 VSUPP IP HD PLUS WHITE wall mount for video camera, white (optional)
- 100 ADAPT BNC RCA adapter (optional)
- 101 ADAPT DVI HDMI adapter (optional)
- 102 ANY WHERE 2 USB over LAN device with 2 USB ports (optional)
- 103 ANY WHERE 5 USB over LAN device with 5 USB ports (optional)
- 104 USB CAM JOYSTICK to control cameras (optional)
- 105 24-key CAM PAD numeric keypad with USB connection for camera control (optional)
- 106 Numeric keypad CAM PAD 48 by 48 keys with USB connection for camera control (optional)
- 107 Connecting cable CB BNBN1 (BNC to BNC) (optional)
- 108 Connecting cable CB BNRC1 (BNC - RCA (Red) (optional)
- 109 Connecting cable CB BNRC5 (BNC - RCA (Black) (optional)
- 110 CB DVI monitor cable (optional)
- 112 Cable HDMI - HDMI MINI for video transmission, length 2m/5m/10m (optional)
- 113 HDMI to HDMI cable for high-definition video transmission, 2m/5m/15m (optional)

- 115 Audio cable CB JACK RCRC with 6.3 mm plugs - double RCA (optional)
- 116 CB MICROPHONE cable (optional)
- 117 CB MON cable for monitor (D-Sub 15 pin - D-Sub 15 pin) (optional)
- 118 CB RCRC JMINI, RCA male to MiniJack 3.5mm 2-wire audio cable (optional)
- 119 Two-wire audio cable CB RCRC PAIR, RCA male to double RCA (as required)
- 120 CB VIDEO video cable (optional)
- 121 CB IP HD PLUS 15 cable for controlling IP cameras (optional)
- 122 CB IP HD PLUS DC 5 cable for connecting IP HD PLUS cameras to each other (optional)
- 123 CB IP HD PLUS EXTENDER adapter for IP HD PLUS signal transmission over Cat5e/6 cable, pair (optional)
- 124 HDMI BOOSTER self-powered HDMI amplifier (optional)
- 125 CARSTAR mobile trolley (on request)
- 126 FIX TOR bracket for attaching the isolation transformer to the trolley (optional)
- 127 Mouse pad small/regular/large (optional)
- 128 Trolley CAR PLUS, mobile, adjustable in height (as required)
- 129 FIX CART plates for fixing trolley legs (optional)
- 130 Power cable kit KIT CB POWER for trolley (optional)
- 131 Set of labels KIT LABEL PLUS for trolley (as required)
- 132 SUPP KEYB PLUS keyboard drawer for trolley (optional)
- 133 SUPP PC PLUS trolley mount (optional)
- 134 SUPP PC PLUS SFF small computer trolley mount (optional)
- 135 Horizontal bar BAR POLE NEW for attaching the bracket to the trolley (optional)
- 136 CAR TORO trolley with shelves, standard size (as required)
- 137 CAR TORO SMALL trolley with shelves, small size (on request)
- 138 COVER TORO back cover for standard CARTORO trolley (optional)
- 139 Rear cover COVER TORO SMALL for small trolley CAR TORO SMALL (optional)
- 140 Isolation transformer TOR 500, Pmax 500 VA, black (on request)
- 141 Isolation transformer TOR 500 W, Pmax 500 VA, white (on request)
- 142 Isolation transformer TR 50, 100 VA (as required)
- 143 NET ISOLAT network isolator for 10/100 Mbit networks, with RJ45 connectors (optional)
- 144 Network cable NET RJ45 with RJ45 connector, length 3m/5m/10m (optional)
- 145 NET SWITCH 4 4-port network switch (optional)
- 146 NET SWITCH 8 8-port network switch (optional)
- 147 NET SWITCH PRO professional network switch (optional)

- 148 Standard Bluetooth USB interface (optional)
- 149 BS USB READER Compact Flash card reader (optional)
- 150 Extension LTM BUTTON EXT for patient button LTM BUTTON TP, length 120 cm (optional)
- LTM BUTTON TP Patient Button for use with SD LTM Series EEG Head Units for Long-Term Monitoring
- 152 PLIERS TP forceps for removing electrodes from the electrode box (optional)
- 153 USB extension cable USB EXT CB (optional)
- 154 USB EXTENDER device for transmitting USB signal over Cat5e/6 cable, including receiver and transmitter (optional)
- 155 BQ CONTACT contact closure control module (optional)
- 156 Patient button with 5 m cable and RCA connector (optional)
- 157 BUTTON BOX for controlling patient buttons (optional)
- 158 Cables for BUTTON BOX (optional)
- 159 Illuminated patient button with 2 m cable and LEMO connector (optional)
- 160 Power supply for BUTTON BOX (optional)
- 161 Cable for oximetry sensor connection unit. Length 3 m (as required)
- 162 Oximetry Sensor Connection Box for use with SD LTM Long-Term EEG Head Units (optional)
- 163 Oximetry sensor "EarClip", adult (optional)
- 164 FingerClip oximetry sensor, adult (optional)
- 165 FingerClip oximetry sensor, child (optional)
- 166 Fixing tape for SD OXY RA, adult size (on request)
- 167 Fixing tape SD OXY RC, child size (as required)
- 168 Flexible oximetry probe, adult size (optional)
- 169 Flexible oximetry sensor, child size (optional)
- 170 Large reusable soft oximeter probe (digital readout 12.5 - 25.5 mm) (optional)
- 171 Medium size reusable soft oximeter probe (digital readout 10 - 19 mm) (optional)
- 172 Reusable soft oximeter probe small size (digital display 7.5 - 12.5 mm) (optional)
- 173 COMP DVI RX KVM extender compact, for DVI + USB HID reception (optional)
- 174 KVM Extender COMP DVI RX AUD compact, for receiving DVI + USB HID + analog audio signal (optional)
- 175 KVM Extender COMP DVI TX compact, for DVI + USB HID transmission (optional)
- 176 COMP DVI TX AUD Compact KVM Extender for DVI + USB HID + analog audio transmission (optional)
- 177 KVM DVI module RX for DVI + USB reception HID (on demand)
- 178 KVM DVI module RX AUD dual for DVI + USB reception HID + analog audio signal (optional)

- 179 KVM DVI module TX for DVI + USB transmission HID (on demand)
- 180 KVM DVI module TX AUD dual for DVI + USB transmission HID + analog audio signal (optional)
- 181 KVM DVI module TX AUD DCAT dual for DVI + USB transmission HID + analog audio signal with two-channel transmission (as required)
- 182 KVM DVI module TX DCAT to DVI + USB transmission HID ,dual transmission (optional)
- 183 KIT KVM KVM module including transmitter and receiver for keyboard, mouse and monitor. (optional)
- 184 KVM module KIT KVM DUAL, including transmitter and receiver, for keyboard, mouse and two monitors (optional)
- 185 KVM - KVM Chassis BODY 2 for 2 DVI modules (optional)
- 186 KVM - KVM Chassis BODY 4 for 4 DVI modules (optional)
- 187 KVM - KVM Chassis BODY 4 R for 4 DVI modules , including rack mounting bracket (optional)
- 188 KVM - KVM Chassis BODY 6 for 6 DVI modules (optional)
- 189 KVM - KVM Chassis BODY 6 R for 6 DVI modules , including rack mounting bracket (optional)
- 190 IP KVM Module, KVM IP RX. Receiver (optional)
- 191 KVM IP, KVM IP TX, Transmitter (optional)
- 192 KVM matrix KVM SWITCH for 8/16/32/48/64/80 ports (on demand)
- 193 Additional SUPP TOPHAT mount for JB LTM TOPHAT Electrode box(optional) 194 HBC EXTERNAL BATTERY cable for external battery connection (optional)
194. Electrode cap (as needed):
- MRINEO 19 X
 - MRINEO 31 X
 - MRINEO 61 X
195. Cable for EL CFM electrodes (as required):
- CB CFM050
 - CB CFM150
196. Electrodes with cable for LTM (as required):
- CB EL TPCO 19
 - CB EL TPCO 31
197. Electrode (plate) EL CFM for 2-channel EEG (on demand)
198. Disk electrode with cable (as required):
- EL TPCO
 - EL TPCO 6
 - EL TPCO G06
 - EL TPCO G10
 - EL TPCO S 100
 - EL TPCO S 150
 - EL TPCO S 200
 - EL TPCO S 250
 - EL TPCO SG 100
 - EL TPCO SG 150

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200. Ear electrode (clip) EL TPEAR with cable (on request)

201. EL TPLINK cable for combining two ear electrodes EL TPEAR (on request)

202. Connecting cable for LTM / LTM 256 (on request)

203. CB PUD150 cable for EL PUD electrodes (optional)

204. Disposable electrodes EL PUD (on request)

205. Electrode monopolar needle (on demand):

- EMG MBOT2530
- EMG MBOT2540
- EMG MBOT3745
- EMG MBOT5050
- EMG MBOT7570

206. Snoring sensor (optional):

- EPM 1200
- EPM 1250

207. EPM 1312 breathing sensor kit, includes sensor (EPM 1387) and 4 tapes in different sizes (EPM 1340 + EPM 1341 + EPM 1342 + EPM 1343) (optional)

208. Pediatric breathing sensor EPM 1390 (on request)

209. Pediatric airflow sensor EPM 1400 (optional)

210. EPM 1439FM/ GM / HM disposable breath sensor kit, includes: EPM 1415 cable and 10 sensors adult 1426 EPMs / medium size 1427 EPMs / infant size 1428 EPMs) (as required)

211. Air flow sensor EPM 1450 (optional)

212. Set of air flow and pressure sensors EPM 1475 (optional)

213. Air pressure and snoring sensor EPM 14830 (optional)

214. Body position sensor EPM 1563 (optional)

215. EPM 1715 motion sensor kit (including piezo transducer (EPM 1770) and pairs of ankle straps (EPMS 1730) (optional)

216. Chest force sensor kit, including interface cable and 4 chest straps (as required):

- EPM 9102 A for adults
- EPM 9102 C pediatric

217. Abdominal force sensor kit, including interface cable and 4 abdominal bands (as required):

- EPM 9101 A for adults
- EPM 9101 C pediatric

218. A set of abdominal and thoracic force sensors, including abdominal and thoracic interface cables, and 4 tapes (as required):

- EPM 914x A for adults
- EPM 914x C pediatric
- EPM 915x A for adults
- EPM 915x C pediatric
- 219. Tape for the breath sensor (as needed):
 - EPMs 1340
 - EPMs 1341
 - EPMs 1342
 - EPMs 1343
- 220. Piezocrystalline breath sensor EPMs 1387 (optional)
- 221. Cable for single-channel airflow sensor EPMs 1415 (optional)
- 222. Thermal breath sensor (on demand):
 - EPMs 1426 adult size
 - EPMs 1427 medium size
 - EPMs 1428 pediatric
- 223. EPMS 1450-B interface for EPM 1450 airflow sensor (optional)
- 224. EPMS 1450- S transducer for EPM 1450 air flow sensor (optional)
- 225. Combined thermistor pressure and airflow sensor EPMs 1475
- 226. CO 2 cannula (as required):
 - EPMs 14802, nasal and oral, male adult
 - EPMs 14803, male for newborns
 - EPMs 14804, male pediatric
 - EPMs 14805, male adult
- 227. Piezoelectric pressure transducer for air flow and snoring EPMs 14830 (on request)
- 228. Elastic tape EPMs 1550 for body position sensor EPM 1520 (optional)
- 229. Converter EPMs 1563 for body position sensor EPM 1563 (optional)
- 230. ThermoCan cannula (optional):
 - EPMs 15802
 - EPMs 15807
- 231. Pair of EPMs 1730 ankle straps (optional)
- 232. Piezoelectric transducer EPMs 1770 for motion sensor EPM 1715 (optional)
- 233. A pair of induction tapes (as needed):
 - 9002-L120 EPMs
 - 9002-L133 EPMs
 - 9002-L40 EPMs
 - 9002-L60 EPMs
 - 9002-L90 EPMs
- 234. Adjustable induction tape EPMs 9010 (on demand)
- 235. Induction interface cable (as needed):
 - EPMs 9101 for the abdomen
 - EPMs 9102 for the thoracic region
- 236. Disposable ECG electrode ECG EL (on demand)
- 237. Cable for ECG electrode (optional):
 - ECG EL TPCLB black
 - ECG EL TPCLR red

IV. ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL STUDIES MICROMED , IN PERFORMANCE: SD LTM 64 PLUS , CONSISTING OF

System software PLUS EVOLUTION - BRAIN QUICK including the following software options:

- 2.1 Software option for recording and viewing EEG, 2-channel CFM version (on request)
- 2.2 Software option for recording and viewing EEG, 8-channel CFM version (on request)
- 2.3 Software option for recording and viewing EEG, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
- 2.4. Software option for recording and viewing digital video-EEG up to 128 channels, version: CLINI C (on request)
- 2.5. Software option for recording and viewing digital video-EEG up to 256 channels, CLINI C version (on request)
- 2.6 Software option for recording and viewing digital video-EEG up to 64 channels, version CLINI C (on request)
- 2.7 Software option for recording and viewing digital video-EEG, version CLINI C (on request)
- 2.8. Software option for recording and viewing digital video-EEG in Full HD format (on request)
- 2.9 Software option for recording and viewing digital video-EEG over the network, version CLINI C (on request)
- 2.10 Software option for recording and viewing digital video-EEG over a network of up to 128 channels, version CLINI C (on request)
- 2.11 Software option for recording and viewing digital video-EEG over a network up to 256 channels, CLINI C version (on request)
- 2.12 Software option for recording and viewing digital video-EEG over a network of up to 64 channels, version CLINI C (on request)
- 2.13 Software option for real-time logging and viewing over the network (as required)
- 2.14 Software option for EEG analysis in "on-line" and "off-line" modes, version: full / CFM / DSA. (on demand)
- 2.15 Software option for transmitting EEG data via TCP/IP (optional)
- 2.16 Software option "ERP -manager" (optional)
- 2.17 Software option "Management of outpatient registrar" (as required)
3. Software options for viewing:
- 3.1. EEG software viewing option, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
- 3.2 Software option for viewing EEG, version: EASY (on request)
- 3.3 Software option for viewing EEG in "off-line" mode, version: full, CFM, DSA (on request)
- 3.4 Software option for connecting to databases of external applications (optional)
- 3.5 EEG software option "Event identification" (as required)
- 3.6 EEG software option "Analysis of high-frequency oscillations" (on request)
- 3.7. Software option EEG "Brain Mapping", version: CLINIC, STUDIO (on request)
- 3.8. EEG software option "Data exchange with MRI" (on request)
- 3.9. Software EEG option "Interface HOLBERG SCORE" (optional)
- 3.10. Software EEG option "Audio playback" (optional)
- 3.11. Software EEG option "Averaging in "off-line" mode" (on demand)
- 3.12. Software option for viewing digital video-EEG, version: CLINIC, STUDIO (on request)
- 3.13. Software option for viewing digital video-EEG in Full HD format (on request)
- 3.14. Software option for viewing digital video-EEG in multichannel mode, version: CLINIC (on request)
- 3.15. Software option for advanced management of the archive server (optional)
- 3.16. Program statistics option (optional)
- 3.17. Software option for sleep program interface (optional)
- 3.18. Software option for background data transfer (optional)
- 3.19. Software option for exporting EEG data to the Viewer (optional)

3.20. Software option to enable HIS interface (optional)

4. Software options for sleep analysis:

- 4.1 Software option for analyzing imported curves in EDF format (optional)
- 4.2 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring 64 channels (on request)
- 4.3 Software option for the analysis of curves recorded on the EEG Head Units of the SD PLUS series (on request)

5. Program control options.

- 5.1 Software option "Centralized control of video cameras" on a USB -key (on request)
- 5.2. Software option "Event management" on the USB -key (on demand)

6. USB key for registration and viewing programs (optional)

7. USB dongle for viewers for 1/5/10/25 users (optional)

8. USB dongle for sleep analysis programs (optional)

9 EXTERNAL BATTERY (optional)

10 EXTERNAL CHARGER for external battery (optional)

11 LCD: 19"/22"/24"/27" (optional)

12. Electrode box JB DSUB PLUS with Dual 37-Pin D - Sub Connector to EEG Head Units for SD Series Long-Term Monitoring LTM PLUS

13 JB TOPHAT PLUS electrode box for all channels with dual 37-pin D-Sub connector to EEG headunits for long-term monitoring SD series LTM PLUS

14 EEG head unit for long-term monitoring SD LTM 64 PLUS (64 channels)

15 Electrode box JB LTM 64 PLUS to the EEG head unit for long-term monitoring SD LTM 64 PLUS (64 channels)

16 HBC EXTERNAL BATTERY cable for connecting an external battery (optional)

17 SD LTM STIM Cortical Stimulator to EEG Head Units for Long-Term Monitoring SD Series LTM (64 channels)

18 LED photostimulator with bracket

19 BOX TERMINAL MULTI interface module for connecting the EEG head unit for long-term monitoring SDLTMPLUS (up to 256 channels)

20 BQ NET POE PLUS interface module for connecting EEG head unit for long-term monitoring SD LTM PLUS , includes SUPP wall mount BQ NET POE 2

21 BQ NET POWER POE power supply for BQ NET POE series interface modules (optional)

22 BQ power supply POWER MULTI for BQ interface module USB MULTI

BQ interface module USB MULTI to EEG headunit for long-term monitoring SD LTM PLUS (up to 256 channels)

BQ interface module USB PLUS to EEG headunit for long-term SD monitoring LTM PLUS (64 channels)

25 USB adapter to serial port ADAPT USB SER (optional)

- 26 Isolated USB to Serial Adapter ADAPT USB SER (optional) ISO
- 27 CB LTM cable for connecting SD cortical stimulator LTM STIM to EEG head unit for long-term SD monitoring LTM64 _
- 28 CB LTM DC cable for connecting DC channels to EEG Head Units for long-term SD monitoring LTM (on demand)
- 29 CB MAG RAPID cable for connecting the Micromed system to Magstim magnetic stimulation devices (D-Sub 26 pin (optional)
- 30 Trigger cable CB RCRC4N for stimulator (optional)
- 31 Trigger cable CB RCRC4R for stimulator (optional)
- 32 Trigger cable CB RCRC1R for stimulator (optional)
- 33 Pair of CB RCRC5 PAIR trigger cables for stimulator (optional)
- 34 CB SER RCA EXT cable for synchronization with external stimulators
- 35 Fiber optic cable 10m/15m/20m/30m/40m/50m/100m (optional)
- 36 HBC EXPRESS cable for connecting an interface board or an interface module to the EEG Head Units (as required)
- 37 HBC INT LTM R cable for connecting an interface board or interface module to EEG Head Units for long-term monitoring SD series LTM
- 38 HBC LTM PLUS cable to EEG headunit for long-term monitoring SD LTM PLUS in 128-channel configuration
- 39 HBC MULTI cable to EEG headunit for long-term monitoring SD LTM PLUS in 256-channel configuration
- 42 LED POWER power supplies for LED photostimulator
- 43 Extension cable LED EXT for LED photostimulator (on request)
- 44 TERM LTM Quick Disconnect Cable to EEG Head Units (optional)
- 45 TERM PLUS quick-disconnect cable to EEG headunit for long-term monitoring SD LTM PLUS in configuration up to 64 channels (on request)
- 46 Quick disconnect cable TERM MULTI 128 to two EEG Head Units for long-term monitoring SD LTM PLUS in a configuration for 128 channels (on request)
- 47 Quick disconnect cable TERM MULTI 256 to four EEG Head Units for long-term monitoring SD LTM PLUS in a configuration for 256 channels (on request)
- 48 CAMERA IP HD PRO (optional)
- 49 CAMERA HD Full HD (1920 x 1080) with PTZ control and automatic IR switching (on demand) (on demand)
- 50 CAMERA IP HD OUTDOOR (optional)
- 51 CAMERA IP HD PLUS video camera (optional)
- 52 CAMERA LAPTOP with USB 2.0, including microphone (optional)
- 53 CAMERA IP COMPACT video camera (optional)
- 54 GC MM video card with dual monitor outputs (optional)

- 55 HDMI EXTENDER device for transmitting HDMI signal over Cat 5 e /6 cable, including receiver and transmitter (optional)
- 56 IP MIC for IP cameras (optional)
- 57 IP POWER power adapter for IP cameras (optional)
- 58 VLAMP IR auto-on (on demand)
- 59 VMIXER PRO 2-input composite switcher (optional)
- 60 VMIXER PRO professional switcher with 5 HDMI inputs (optional)
- 61 Bracket ARM GM LTM for EEG head units for long-term monitoring of the SD LTM series, complete (on request)
- 62 Vertical bracket ARM LCD DOUBLE for two LCD displays for fixing on the cart (optional)
- 63 Horizontal bracket ARM LCD HORIZ for two LCD displays for fixing on a trolley (optional)
- 64 ARM LCD SW Vertical LCD Bracket for Trolley Fixing (Optional)
- 65 Adjustable BAG BELT HP for fixing the headunit bag to the patient
- 66 Bag BAG LTM 128 for two EEG Head Units for long-term monitoring of the SD LTM series
- 67 BAG LTM HP Bag for EEG Head Unit for Long-Term Monitoring SD LTM Series (optional)
- 68 BASE STAND (optional)
- 69 BQ CLAMP (optional)
- 70 CAM JOIN swivel mount for video camera on a stand (optional)
- 71 CAM STAND for mounting a video camera or head unit on a trolley (optional)
- 72 Fixing FIX GM for trolley bracket (optional)
- 73 Fixture FIX GM TW for table bracket (optional)
- 74 JOIN GM connection for attaching the head unit holder to the bracket (optional)
- 75 LTM SCREW screw for attaching the electrode box to the head unit (optional)
- 76 LED STAND PHOTIC wheel base for LED photostimulator (on request)
- 77 LED SUPP PHOTIC bracket for fixing the LED photostimulator (optional)
- 78 STAND CAM SW video camera holder assembly (optional)
- 79 STAND POLE NEW on wheels with weight (optional)
- 80 SUPP BQ NET POE 2 wall mount for BQ interface modules NET POE2 / BQNETPOEM
- 81 SUPP LTM Holder for SD LTM Series EEG Head Units (optional)
- 82 Holder SUPP LTM 128 for two EEG Head Units for long-term monitoring of the SD LTM series in a configuration for 128 channels (on request)
- 83 Holder SUPP LTM 256 for four EEG Head Units for long-term monitoring of the SD LTM series in a configuration for 256 channels (on request)
- 85 VCLAMP mount for video camera or IR light (optional)
- 86 VSTAND tripod stand for video camera or IR light (optional)

- 88 VSUPP IP metal support for attaching the IP video camera to the swivel mount (optional)
- 89 VSUPP IP IR metal support for attaching the IP camera to the rack (optional)
- 90 VSUPP IP HD PLUS BLACK wall mount for video camera, black (optional)
- 91 VSUPP IP HD PLUS WHITE wall mount for video camera, white (optional)
- 92 ADAPT BNC RCA adapter (optional)
- 93 ADAPT DVI HDMI adapter (optional)
- 94 ANY WHERE 2 USB over LAN device with 2 USB ports (optional)
- 95 ANY WHERE 5 USB over LAN device with 5 USB ports (optional)
- 96 USB CAM JOYSTICK to control cameras (optional)
- 97 24-key CAM PAD numeric keypad with USB connection for camera control (optional)
- 98 Numeric keypad CAM PAD 48 by 48 keys with USB connection for camera control (optional)
- 99 Connecting cable CB BBNB1 (BNC to BNC) (optional)
- 100 Connecting cable CB BNRC1 (BNC - RCA (Red) (as required)
- 101 Connecting cable CB BNRC5 (BNC - RCA (Black) (optional)
- 102 CB DVI monitor cable (optional)
- 103 Cable HDMI - HDMI MINI for video transmission, length 2m/5m/10m (optional) (optional)
- 104 HDMI to HDMI cable for high-definition video transmission, 2m/5m/15m (optional)
- 105 Audio cable CB JACK JMINI with 6.3 mm plugs - MiniJack 3.5 mm (optional)
- 106 Audio cable CB JACK RCRC with 6.3 mm plugs - double RCA (optional)
- 107 CB MICROPHONE cable (optional)
- 108 CB MON cable for monitor (D-Sub 15 pin - D-Sub 15 pin) (optional)
- 109 CB RCRC JMINI, RCA male to MiniJack 3.5mm 2-wire audio cable (optional)
- 110 Two-wire audio cable CB RCRC PAIR, RCA male to double RCA (as required)
- 111 CB VIDEO video cable (optional)
- 112 CB IP HD PLUS 15 cable for controlling IP cameras (optional)
- 113 CB IP HD PLUS DC 5 cable for connecting IP HD PLUS cameras to each other (optional)
- 114 CB IP HD PLUS EXTENDER adapter for IP HD PLUS signal transmission over Cat5e/6 cable, pair (optional)
- 115 HDMI BOOSTER self-powered HDMI amplifier (optional)
- 116 CARSTAR mobile trolley (on request)
- 117 FIX TOR bracket for attaching the isolation transformer to the trolley (optional)
- 118 Mouse pad small/regular/large (optional)
- 119 CAR PLUS trolley, mobile, adjustable in height (as required)
- 120 FIX CART plates for fixing trolley legs (optional)

- 121 Power cable kit KIT CB POWER for trolley (optional)
- 122 KIT LABEL PLUS for trolley (as required)
- 123 SUPP KEYB PLUS keyboard drawer for trolley (optional)
- 124 SUPP PC PLUS trolley mount (optional)
- 125 Small computer mount SUPP PC PLUS SFF on cart (optional)
- 126 Horizontal bar BAR POLE NEW for attaching the bracket to the trolley (optional)
- 127 CAR TORO trolley with shelves, standard size (as required)
- 128 CAR TORO SMALL trolley with shelves, small size (on request)
- 129 COVER TORO rear cover for standard CARTORO trolley (optional)
- 130 Rear cover COVER TORO SMALL for small cart CAR TORO SMALL (optional)
- 131 Isolation transformer TOR 500, Pmax 500 VA, black (on request)
- 132 Isolation transformer TOR 500 W, Pmax 500 VA, white (on request)
- 133 Isolating transformer TR 50, 100 VA (on request)
- 134 NET ISOLAT network isolator for 10/100 Mbit networks, with RJ45 connectors (optional)
- 135 Network cable NET RJ45 with RJ45 connector, length 3m/5m/10m (as required)
- 136 NET SWITCH 4 4-port network switch (optional)
- 137 NET SWITCH 8 8-port network switch (optional)
- 138 NET SWITCH PRO professional network switch (optional)
- 139 Standard Bluetooth USB interface (optional)
- 140 BS USB READER Compact Flash card reader (optional)
- 141 Extension LTM BUTTON EXT for patient button LTM BUTTON TP, length 120 cm (optional)
- LTM BUTTON TP Patient Button for use with SD LTM Series EEG Head Units for Long-Term Monitoring
- 143 PLIERS TP forceps for removing electrodes from the electrode box (optional)
- 144 USB extension cable USB EXT CB (optional)
- 145 USB EXTENDER device for transmitting USB signal over Cat5e/6 cable, including receiver and transmitter (optional)
- 146 BQ CONTACT contact closure control module (optional)
- 147 Patient button with 5 m cable and RCA connector (optional)
- 148 BUTTON BOX for controlling patient buttons (optional)
- 149 Cables for BUTTON BOX (optional)
- 150 Illuminated patient button with 2 m cable and LEMO connector (optional)
- 151 Power supply for BUTTON BOX (optional)
- 152 Cable for connection unit for oximetry sensors. Length 3 m (as required)

- 153 Oximetry Sensor Connection Box for use with EEG Head Units for Long-Term Monitoring SD LTM Series (as required)
- 154 Oximetry sensor "EarClip", adult (optional)
- 155 FingerClip oximetry sensor, adult (optional)
- 156 FingerClip oximetry sensor, child (optional)
- 157 Fixing tape for SD OXY RA, adult size (on request)
- 158 Fixing tape SD OXY RC, child size (as required)
- 159 Flexible oximetry probe, adult size (optional)
- 160 Flexible oximetry sensor, child size (optional)
- 161 Large reusable soft oximeter probe (digital readout 12.5 - 25.5 mm) (optional)
- 162 Reusable medium size soft oximeter probe (digital readout 10 - 19 mm) (optional)
- 163 Reusable soft oximeter probe small size (digital readout 7.5 - 12.5 mm) (optional)
- 164 COMP DVI RX KVM extender compact, for DVI + USB HDI reception (optional)
- 165 KVM Extender COMP DVI RX AUD compact, for receiving DVI + USB HDI + analog audio signal (optional)
- 166 KVM Extender COMP DVI TX compact, for DVI + USB HDI transmission (optional)
- 167 KVM Extender COMP DVI TX AUD compact, for DVI + USB HDI + analog audio transmission (optional)
- 168 KVM DVI module RX for DVI + USB reception HID (on demand)
- 169 KVM DVI module RX AUD dual for DVI + USB reception HID + analog audio signal (optional)
- 170 KVM DVI module TX for DVI + USB transmission HID (on demand)
- 171 KVM DVI module TX AUD dual for DVI + USB transmission HID + analog audio signal (optional)
- 172 KVM DVI module TX AUD DCAT dual for DVI + USB transmission HID + analog audio signal with two-channel transmission (as required)
- 173 KVM DVI module TX DCAT to DVI + USB transmission HID ,dual transmission (optional)
- 174 KIT KVM module including transmitter and receiver for keyboard, mouse and monitor. (optional)
- 175 KVM module KIT KVM DUAL, including transmitter and receiver, for keyboard, mouse and two monitors (optional)
- 176 KVM - KVM Chassis BODY 2 for 2 DVI modules (optional)
- 177 KVM -KVM Chassis BODY 4 for 4 DVI modules (optional)
- 178 KVM - KVM Chassis BODY 4 R for 4 DVI modules , including rack mounting bracket (optional)
- 179 KVM - KVM Chassis BODY 6 for 6 DVI modules (optional)
- 180 KVM -KVM Chassis BODY 6 R for 6 DVI modules , including rack mounting bracket

181M IP module, KVM IP RX, receiver (optional)

182M-module IP, KVM IP TX, transmitter (optional)

183.KVM matrix KVM SWITCH for 8/16/32/48/64/80 ports (on demand)

184. Electrode cap (as needed):

- MRINEO 19 X

- MRINEO 31 X

- MRINEO 61 X

185. Cable for EL CFM electrodes (as required):

- CB CFM050

- CB CFM150

186. Electrodes with cable for LTM (as required):

- CB EL TPCO 19

- CB EL TPCO 31

187. Electrode (plate) EL CFM for 2-channel EEG (on demand)

188. Disk electrode with cable (as required):

- EL TPCO

- EL TPCO 6

- EL TPCO G06

- EL TPCO G10

- EL TPCO S 100

- EL TPCO S 150

- EL TPCO S 200

- EL TPCO S 250

- EL TPCO SG 100

- EL TPCO SG 150

- EL TPCO SG 200

- EL TPCO SG 250

- EL TPCO SLT 100

- EL TPCO SLT 150

- EL TPCO SLT 200

- EL TPCO SLT 250

- EL TPCO SLTG 100

- EL TPCO SLTG 150

- EL TPCO SLTG 200

- EL _ TPCO SLTG 250

189. Pediatric disk electrode with cable (on demand):

- EL TPCO S6 100

- EL TPCO S6 150

- EL TPCO S6 200

- EL TPCO S6 250

- EL TPCO SG6 100

- EL TPCO SG6 150

- EL TPCO SG6 200

- EL TPCO SG6 250

190. Ear electrode (clip) EL TPEAR with cable (on request)

191. EL TPLINK cable for combining two ear electrodes EL TPEAR (optional)

192.Connecting cable for LTM / LTM 256 (on request)

193.CB PUD150 cable for EL PUD electrodes (optional)

194. Disposable electrodes EL PUD (on demand)

195. Electrode monopolar needle (on demand):

- EMG MBOT2540
- EMG MBOT3745
- EMG MBOT5050
- EMG MBOT7570
- 196. Snoring sensor (on demand):
 - EPM 1200
 - EPM 1250
- 197. EPM 1312 breathing sensor kit, includes sensor (EPM 1387) and 4 tapes in different sizes (EPM 1340 + EPM 1341 + EPM 1342 + EPM 1343) (optional)
- 198. Pediatric breathing sensor EPM 1390 (on request)
- 199. Pediatric air flow sensor EPM 1400 (on demand)
- 200. EPM 1439FM/ GM / HM disposable breath sensor kit , includes: EPM 1415 cable and 10 sensors adult 1426 EPMs / medium size 1427 EPMs / infant size 1428 EPMs) (as required)
- 201. Air flow sensor EPM 1450 (optional)
- 202. Set of air flow and pressure sensors EPM 1475 (optional)
- 203. Air pressure and snoring sensor EPM 14830 (optional)
- 204. Body position sensor EPM 1563 (optional)
- 205. EPM 1715 motion sensor kit (including piezoelectric transducer (EPM 1770) and pairs of ankle straps (EPMS 1730) (optional)
- 206. Chest force sensor kit, including interface cable and 4 chest straps (as required):
 - EPM 9102 A for adults
 - EPM 9102 C pediatric
- 207. Abdominal force sensor kit, including interface cable and 4 abdominal bands (as required):
 - EPM 9101 A for adults
 - EPM 9101 C pediatric
- 208. A set of abdominal and thoracic force sensors, including abdominal and thoracic interface cables, and 4 tapes (as required):
 - EPM 914x A for adults
 - EPM 914x C pediatric
 - EPM 915x A for adults
 - EPM 915x C pediatric
- 209. Tape for the breath sensor (as needed):
 - EPMs 1340
 - EPMs 1341
 - EPMs 1342
 - EPMs 1343
- 210. Piezocrystalline breath sensor EPMs 1387 (optional)
- 211. Cable for single-channel airflow sensor EPMs 1415 (optional)
- 212. Thermal breath sensor (as required):
 - EPMs 1426 adult size
 - EPMs 1427 medium size
 - EPMs 1428 pediatric
- 213. EPMS 1450-B interface for EPM 1450 airflow sensor (optional)
- 214. EPMS 1450- S transducer for EPM 1450 airflow sensor (optional)
- 215. Combined thermistor pressure and airflow sensor EPMs 1475
- 216. CO 2 cannula (as required):
 - EPMs 14802, nasal and oral, male adult
 - EPMs 14803, male for newborns
 - EPMs 14804, male pediatric
 - EPMs 14805, male adult
- 217. Piezoelectric pressure transducer for air flow and snoring EPMs 14830 (on request)
- 218. Elastic tape EPMs 1550 for body position sensor EPM 1520 (optional)
- 219. Converter EPMs 1563 for body position sensor EPM 1563 (optional)
- 220. ThermoCan cannula (optional):
 - EPMs 15802
 - EPMs 15807

- 43
- 221. Pair of EPMs 1730 ankle straps (optional)
 - 222. Piezoelectric transducer EPMs 1770 for motion sensor EPM 1715 (optional)
 - 223. A pair of induction tapes (as needed):
 - 9002-L120 EPMs
 - 9002-L133 EPMs
 - 9002-L40 EPMs
 - 9002-L60 EPMs
 - 9002-L90 EPMs
 - 224. Adjustable induction tape EPMs 9010 (on request)
 - 225. Induction interface cable (as required):
 - EPMs 9101 for the abdomen
 - EPMs 9102 for the thoracic region
 - 226. ECG electrode disposable ECG EL (on demand)
 - 227. Cable for ECG electrode (optional):
 - ECG EL TPCLB black
 - ECG EL TPCLR red

V. ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL STUDIES MICROMED , PERFORMED BY: SD PLUS FLEXI CLINIC , CONSISTING OF:

System software PLUS EVOLUTION - BRAIN QUICK including the following software options:

Software option for registration and viewing:

- 2.1 Software option for recording and viewing EEG, 2-channel CFM version (on request)
- 2.2 Software option for recording and viewing EEG, 8-channel CFM version (on request)
- 2.3 Software option for recording and viewing EEG, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
- 2.4 Software option for recording and viewing digital video-EEG, version CLINIC (on request)
- 2.5. Software option for recording and viewing digital video-EEG in FullHD format (on request)
- 2.6 Software option for recording and viewing digital video-EEG over the network, version CLINIC (on request)
- 2.7 Software option for real-time logging and viewing over the network (as required)
- 2.8 Software option for EEG analysis in "on-line" and "off-line" modes, version: full / CFM / DSA. (on demand)
- 2.9 Software option for transmitting EEG data via TCP/IP (optional)
- 2.10 Software option " ERP -manager" (optional)
- 2.11 Software option "Management of outpatient registrar" (as required)
- 3. Software options for viewing:
 - 3.1. EEG software viewing option, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
 - 3.2 Software option for viewing EEG, version: EASY (on request)
 - 3.3 Software option for viewing EEG in "off-line" mode, version: full, CFM, DSA (on request)
 - 3.4 Software option for connecting to databases of external applications (optional)
 - 3.5 EEG software option "Event identification" (as required)
 - 3.6 EEG software option "Analysis of high-frequency oscillations" (on request)

- 3.7. Software option EEG "Brain Mapping", version: CLINIC, STUDIO (on request)
- 3.8. EEG software option "Data exchange with MRI" (on request)
- 3.9. Software EEG option "Interface HOLBERG SCORE" (optional)
- 3.10. Software EEG option "Audio playback" (optional)
- 3.11. Software EEG option "Averaging in "off-line" mode" (on demand)
- 3.12. Software option for viewing digital video-EEG, version: CLINIC, STUDIO (on request)
- 3.13. Software option for viewing digital video-EEG in Full HD format (on request)
- 3.14. Software option for viewing digital video-EEG in multichannel mode, version: CLINIC (on request)
- 3.15. Software option for advanced management of the archive server (optional)
- 3.16. Program statistics option (optional)
- 3.17. Software option for sleep program interface (optional)
- 3.18. Software option for background data transfer (optional)
- 3.19. Software option for exporting EEG data to the Viewer (optional)
- 3.20. Software option to enable HIS interface (optional)
4. Software options for sleep analysis:
- 4.1 Software option for analyzing imported curves in EDF format (optional)
- 4.2 Software option for the analysis of curves recorded on the EEG Head Units of the SD PLUS series (on request)
5. Program control options.
- 5.1 Software option "Centralized control of video cameras" on a USB -key (on request)
- a. Software option "Event management" on a USB key (optional)
6. USB key for registration and viewing programs (optional)
7. USB dongle for viewers for 1/5/10/25 users (optional)
8. USB dongle for sleep analysis programs (optional)
- BQUSBPLUS interface module to EEG head unit for long-term monitoring SDLTMPPLUS (64 channels)
- 10 IIBC EXPRESS cable for connecting an interface board or interface module to EEG Head Units (as required)
- 11 LCD: 19"/22"/24"/27" (optional)
- 12 TERM EXPRESS quick-disconnect cable to the EEG Head Unit for long-term monitoring SD LTM EXPRESS in a configuration of up to 64 channels (on request)
- 13 EEG head unit SD PLUS CLINIC (29 channels)
- 14 LED photostimulator with bracket
- 15 BQ NET POE 2 interface module for connecting SD PLUS EEG Headunit, includes SUPP BQ NET POE 2 wall mount and BQ interface module USB EXPRESS
16. HBC BT/RJ45 extension adapter for BOX TERMINAL EXPRESS interface module (module side) (optional)
- 17 BQ NET POWER POE power supply for BQ NET POE series interface modules (optional)
- 18 . BQ PCI EXPR LTM FULL interface board based on PCI Express FULL Size to control SD PLUS EEG headunit
- 19 BQ PCI EXPR LTM HALF interface board based on PCI Express HALF Size to control

- 20 BQ USB EXPRESS Interface Module to SD PLUS EEG Head Unit
- 21 TERMINAL interface module EXPRESS to Head Units EEG SD PLUS
- 22 USB adapter to serial port ADAPT USB SER (optional)
- 23 Isolated USB to Serial Adapter ADAPT USB SER (optional) ISO
- 24 CB MAG RAPID cable for connecting the Micromed system to Magstim magnetic stimulation devices (D-Sub 26 pin (optional)
- 25 Trigger cable CB RCRC4N for stimulator (optional)
- 26 Trigger cable CB RCRC4R for stimulator (optional)
- 27 Trigger cable CB RCRC1R for stimulator (optional)
- 28 Pair of CB RCRC5 PAIR trigger cables for stimulator (optional)
- 29 CB SER RCA EXT cable for synchronization with external stimulators
- 30 Fiber optic cable 10m/15m/20m/30m/40m/50m/100m (optional)
- 31 LED POWER power supply for LED photostimulator
- 32 180 Extension cable LED EXT for LED photostimulator (on request)
- 33 CAMERA IP HD PRO video camera (optional)
- 34 CAMERA HD Full HD (1920 x 1080) camcorder with PTZ control and automatic IR switching (on demand)
- 35 CAMERA IP HD OUTDOOR video camera (optional)
- 36 CAMERA IP HD FullHD (1920 x 1080) camera with PTZ control and automatic IR switching for IP video (on demand)
- 37 Camcorder CAMERA LAPTOP with USB 2.0, including microphone (optional)
- 38 CAMERA IP COMPACT video camera (optional)
- 39 GC MM video card with dual monitor output (optional)
- 40 HDMI EXTENDER device for transmitting HDMI signal over Cat5e/6 cable, including receiver and transmitter (optional)
- 41 IP MIC microphone for IP video cameras (optional)
- 42 IP POWER power adapter for IP video cameras (optional)
- 43 VLAMP IR illuminator with auto-on (on demand)
- 44 Composite switcher VMIXER PRO with 2 inputs (optional)
- 45 Professional switcher VMIXER PRO with 5 HDMI inputs (optional)
- 46 Bracket ARM GM PLUS for Head Units EEG series SD PLUS, assembled (on request)
- 47 Vertical bracket ARM LCD DOUBLE for two LCD displays for fixing on the cart (optional)
- 48 Horizontal bracket ARM LCD HORIZ for two LCD displays for fixing on the cart (optional)
- 49 Vertical bracket ARM LCD SW for LCD displays for fixing on the trolley (optional)
- 50 HBC HT/RJ45 extension adapter for BOX TERMINAL EXPRESS interface module (PC side)(as required)
- 51 Rack with base BASESTAND (optional)
- 52 BQ CLAMP (optional)
- 53 Swivel mount CAM JOIN for a video camera on a stand (optional)
- 54 CAM STAND for mounting a video camera or head unit on a trolley (optional)
- 55 FIX GM mount for trolley bracket (optional)
- 56 Fixing FIX GM TW for table bracket (optional)
- 57 JOIN GM connection for attaching the head unit holder to the bracket (optional)
- 58 SUPP PLUS Holder for SD PLUS EEG Head Units (optional)
- 59 LED STAND PHOTIC wheel base for LED photostimulator (on demand)
- 60 Bracket LED SUPP PHOTIC for fixing the LED photostimulator (optional)
- 61 Holder STAND CAM SW video camera assembly (optional)
- 62 Stand STAND POLE NEW on wheels with weighting (optional)
- 63 SUPP BQ NET POE2 wall mount for BQ NET POE 2 interface modules BQ NET POE 2 / BQ NET POE M (optional) (optional)
- 64 VCLAMP mount for video camera or IR light (optional)
- 65 VSTAND tripod stand for camcorder or IR illuminator (optional)

- 66 VSUPP Wall Mount for Video Camera (Optional)
- 67 VSUPP IP metal support for attaching an IP video camera to a swivel mount (optional)
- 68 VSUPP IP IR metal support for attaching the IP video camera to the rack (optional)
- 69 Wall mount VSUPP IP HD PLUS BLACK for video camera, black (optional)
- 70 VSUPP IP HD PLUS WHITE wall mount for video camera, white (optional)
- 71 ADAPT BNC RCA adapter (optional)
- 72 ADAPT DVI HDMI adapter (optional)
- 73 ANY WHERE 2 USB over LAN device with 2 USB ports (optional)
- 74 ANY WHERE 5 USB over LAN device with 5 USB ports (optional)
- 75 USB joystick CAM JOYSTICK for camera control (optional)
- 76 24-key CAM PAD numeric keypad with USB connection for camera control (optional)
- 77 Numeric keypad CAM PAD 48 x 48 keys with USB connection for camera control (optional)
- 78 Connecting cable CB BNBNI (BNC to BNC) (optional)
- 79 Connecting cable CB BNRC1 (BNC - RCA (Red) (as required)
- 80 Connecting cable CB BNRC5 (BNC - RCA (Black) (optional)
- 81 CB DVI cable for monitor (optional)
- 82 HDMI - HDMI MINI cable for video transmission, length 2m/5m/10m (optional)
- 83 HDMI to HDMI cable for high-definition video transmission, 2m/5m/15m long (optional)
- 84 Audio cable CB JACK JMINI with 6.3 mm plugs - MiniJack 3.5 mm (optional)
- 85 Audio cable CB JACK RCRC with 6.3 mm plugs - double RCA (optional)
- 86 Microphone cable CB MICROPHONE (optional)
- 87 CB MON cable for monitor (D-Sub 15 pin - D-Sub 15 pin) (optional)
- 88 Two-wire audio cable CB RCRC JMINI, RCA male-MiniJack 3.5 mm (optional)
- 89 Two-wire audio cable CB RCRC PAIR, RCA male-double RCA (as required)
- 90 CB VIDEO video cable as required)
- 91 CB IP HD PLUS 15 cable for controlling IP cameras (optional)
- 92 CB IP HD PLUS DC 5 cable for connecting IP HD PLUS cameras to each other (optional)
- 93 CB IP HD PLUS EXTENDER adapter for IP HD PLUS signal transmission over Cat5e/6 cable, pair (optional)
- 94 HDMI BOOSTER self-powered HDMI amplifier (optional)
- 95 Trolley CARSTAR mobile (on demand)
- 96 Bracket FIX TOR for attaching the isolation transformer to the trolley (optional)
- 97 Mouse pad small/regular/large (as needed)
- 98 Trolley CAR PLUS mobile, adjustable in height (as required)
- 99 FIX CART plates for fixing trolley legs (optional)
- 100 Power cable set KIT CB POWER for trolley (as required)
- 101 Set of labels KIT LABEL PLUS for the trolley (as required)
- 102 SUPP KEYB PLUS keyboard drawer for trolley (optional)
- 103 SUPP PC PLUS trolley mount (optional)
- 104 Small computer mount SUPP PC PLUS SFF on cart (optional)
- 105 Horizontal bar BAR POLE NEW for attaching the bracket to the trolley (optional)
- 106 Trolley CAR TORO with shelves, standard size (as required)
- 107 Trolley CAR TORO SMALL with shelves, small size (on demand)
- 108 COVER TORO rear cover for standard CARTORO trolley (optional)
- 109 Back cover COVER TORO SMALL for small cart CAR TORO SMALL (as required)
- 110 Isolation transformer TOR 500, Pmax 500 VA, black (on request)
- 111 Isolation transformer TOR 500 W, Pmax 500 VA, white (on request)
- 112 Isolation transformer TR 50, 100 VA (on request)
- 113 NET ISOLAT network isolator for 10/100 Mbps networks, with RJ45 connectors (optional)
- 114 Network cable NET RJ45 with RJ45 connector, length 3m/5m/10m (optional)
- 115 4-port network switch NET SWITCH 4 (optional)
- 116 8-port network switch NET SWITCH 8 (optional)

- 117 Professional network switch NET SWITCH PRO (optional)
- 118 Standard Bluetooth USB interface (optional)
- 119 BS USB READER for Compact Flash card reader (optional)
- 120 PLIERS TP forceps for removing electrodes from the electrode box (optional)
- 121 USB extension cable USB EXT CB (optional)
- 122 USB EXTENDER device for transmitting a USB signal over a Cat5e / 6 cable, including a receiver and transmitter (optional)
- 123 BQ CONTACT contact closure control module (optional)
- 124 Patient button with 5m cable and RCA connector (optional)
- 125 BUTTON BOX to control patient buttons (optional)
- 126 Cables for BUTTON BOX (optional)
- 127 Illuminated patient button with 2 m cable and LEMO connector (optional)
- 128 Power supply for BUTTON BOX (optional)
- 129 Cable for connection unit for oximetric sensors. Length 3 m (as required)
- 130 Oximetry Sensor Connection Box for use with SD PLUS Series EEG Head Units
- 131 Oximetric sensor "EarClip". adult (on request)
- 132 FingerClip oximetry sensor. adult (optional)
- 133 Oximetric sensor "FingerClip", child (on request)
- 134 Fixing tape for SD OXY RA, adult size (as required)
- 135 Fixing tape SD OXY RC, child size (as required)
- 136 Flexible oximetry probe. adult size (optional)
- 137 Flexible oximetry sensor. child size (optional)
- 138 Large size reusable soft oximeter probe (digital readout 12.5 - 25.5 mm) (optional)
- 139 Reusable medium size soft oximeter probe (digital readout 10 - 19 mm) (optional)
- 140 Reusable soft oximetry sensor small size (digital display 7.5 - 12.5 mm) (optional)
- 141 KVM-extender COMP DVI RX compact. for receiving DVI + USB HDI (optional)
- 142 KVM Extender COMP DVI RX AUD compact, for receiving DVI + USB HDI + analog audio signal (optional)
- 143 KVM Extender COMP DVI TX compact, for DVI + USB HDI transmission (optional)
- 144 COMP DVI TX AUD compact KVM extender for DVI + USB HDI + analog audio transmission (optional)
- 145 DVI RX KVM module for DVI + USB HID reception (optional)
- 146 KVM Module DVI RX AUD Dual to receive DVI + USB HID + analog audio signal (as required)
- 147 DVI TX KVM module for DVI + USB HID transmission (optional)
- 148 KVM Module DVI TX AUD Dual to Transmit DVI + USB HID + Analog Audio Signal (As Required)
- 149 KVM Module DVI TX AUD DCAT Dual for DVI + USB HID + Analog Audio Signal Transmission, Dual Channel Transmission (optional)
- 150 DVI TX DCAT KVM Module for DVI + USB HID Dual Link Transmission (optional)
- 151 KIT KVM KVM module including transmitter and receiver for keyboard, mouse and monitor. (optional)
- 152 KVM module KIT KVM DUAL, including transmitter and receiver, for keyboard, mouse and two monitors (optional)
- 153 KVM Chassis KVM BODY 2 for 2 DVI modules (optional)
- 154 KVM Chassis KVM BODY 4 for 4 DVI modules (optional)
- 155 KVM BODY 4R Chassis for 4 DVI modules, including rack mount bracket (optional)
- 156 KVM Chassis KVM BODY 6 for 6 DVI modules (optional)
- 157 KVM BODY 6R Chassis for 6 DVI modules, including rack mount bracket (optional)
- 158 KVM IP Module, KVM IP RX. Receiver (optional)
- 159 KVM IP module, KVM IP TX, Transmitter (optional)
- 160 KVM matrix KVM SWITCH for 8/16/32/48/64/80 ports (on demand)

161. Electrode cap (optional):

- MRINEO 19 X
- MRINEO 31 X
- MRINEO 61 X

162. Cable for EL CFM electrodes (as required):

- CB CFM050
- CB CFM150

163. Electrodes with cable for LTM (as required):

- CB EL TPCO 19
- CB EL TPCO 31

164. Electrode (plate) EL CFM for 2-channel EEG (on demand)

165. Disk electrode with cable (as needed):

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167. Ear electrode (clip) EL TPEAR with cable (on request)

168. EL TPLINK cable for combining two ear electrodes EL TPEAR (optional)

169. Connecting cable for LTM / LTM 256 (on request)

170. CB PUD150 cable for EL PUD electrodes (optional)

171. Disposable electrodes EL PUD (on demand)

172. Electrode monopolar needle (on demand):

- EMG MBOT2530
- EMG MBOT2540
- EMG MBOT3745

- EMG MBOT7570
- 173. Snoring sensor (optional):
 - EPM 1200
 - EPM 1250
- 174. EPM 1312 breathing sensor kit, includes sensor (EPM 1387) and 4 tapes in different sizes (EPM 1340 + EPM 1341 + EPM 1342 + EPM 1343) (optional)
- 175. Pediatric breathing sensor EPM 1390 (on demand)
- 176. Pediatric air flow sensor EPM 1400 (on demand)
- GM / HM disposable breath sensor set , includes: EPM 1415 cable and 10 sensors adult 1426 EPMs / medium size 1427 EPMs / infant size 1428 EPMs) (as required)
- 178. Air flow sensor EPM 1450 (optional)
- 179. Set of air flow and pressure sensors EPM 1475 (optional)
- 180. Air pressure sensor and snoring EPM 14830 (optional)
- 181. Body position sensor EPM 1563 (optional)
- 182. EPM 1715 motion sensor kit (includes piezo transducer (EPM 1770) and pairs of ankle straps (EPMS 1730) (optional)
- 183. Chest force sensor kit, including interface cable and 4 chest straps (as required):
 - EPM 9102 A for adults
 - EPM 9102 C pediatric
- 184. Set of abdominal force sensor, including interface cable and 4 abdominal bands (as required):
 - EPM 9101 A for adults
 - EPM 9101 C pediatric
- 185. A set of abdominal and thoracic force sensors, including abdominal and thoracic interface cables, and 4 tapes (as required):
 - EPM 914x A for adults
 - EPM 914x C pediatric
 - EPM 915x A for adults
 - EPM 915x C pediatric
- 186. Tape for the breath sensor (as needed):
 - EPMs 1340
 - EPMs 1341
 - EPMs 1342
 - EPMs 1343
- 187. Piezocrystalline breath sensor EPMs 1387 (optional)
- 188. Cable for single-channel airflow sensor EPMs 1415 (optional)
- 189. Thermal breath sensor (on demand):
 - EPMs 1426 adult size
 - EPMs 1427 medium size
 - EPMs 1428 pediatric
- 190. EPMs 1450-B interface for EPM 1450 airflow sensor (optional)
- 191. EPMS 1450- S transducer for EPM 1450 airflow sensor (optional)
- 192. Combined thermistor pressure and airflow sensor EPMs 1475
- CO 2 cannula (as required):
 - EPMs 14802, nasal and oral, male adult
 - EPMs 14803, male for newborns
 - EPMs 14804, male pediatric
 - EPMs 14805, male adult
- 194. Piezoelectric pressure transducer for airflow and snoring EPMs 14830 (optional)
- 195. Elastic tape EPMs 1550 for body position sensor EPM 1520 (optional)
- 196. EPMs 1563 transducer for EPM 1563 body position sensor (optional)
- 197. ThermoCan cannula (as required):
 - EPMs 15802
 - EPMs 15807
- 198. Pair of EPMs 1730 ankle straps (optional)
- 199. Piezoelectric transducer EPMs 1770 for motion sensor EPM 1715 (optional)

200. A pair of induction tapes (as needed):
- 9002-L120 EPMs
 - 9002-L133 EPMs
 - 9002-L40 EPMs
 - 9002-L60 EPMs
 - 9002-L90 EPMs
201. Adjustable induction tape EPMs 9010 (as required)
202. Induction interface cable (as required):
- EPMs 9101 for the abdomen
 - EPMs 9102 for the thoracic region
203. Disposable ECG electrode ECG EL (optional)
204. Cable for ECG electrode (optional):
- ECG EL TPCLB black
 - ECG EL TPCLR red

VI. ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL STUDIES MICROMED , PERFORMED BY: SD PLUS FLEXI , COMPOSED

System software PLUS EVOLUTION - BRAIN QUICK, including the following software options:

Software option for registration and viewing:

- 2.1 Software option for recording and viewing EEG, 2-channel CFM version (on request)
- 2.2 Software option for recording and viewing EEG, 8-channel CFM version (on request)
- 2.3 Software option for recording and viewing EEG, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
- 2.5 Software option for recording and viewing digital video-EEG up to 64 channels, version CLINI C (on request)
- 2.6 Software option for recording and viewing digital video-EEG, version CLINI C (on request)
- 2.7. Software option for recording and viewing digital video-EEG in Full HD format (on request)
- 2.8 Software option for recording and viewing digital video-EEG over the network, version CLINI C (on request)
- 2.9 Software option for recording and viewing digital video-EEG over a network of up to 64 channels, version CLINI C (on request)
- 2.10 Software option for logging and viewing over the network in real time (as required)
- 2.11 Software option for EEG analysis in "on-line" and "off-line" modes, version: full / CFM / DSA. (on demand)
- 2.12 Software option for transmitting EEG data via TCP/IP (optional)
- 2.13 Software option " ERP -manager" (optional)
- 2.14 Software option "Management of the outpatient registrar" (as required)
3. Software options for viewing:
 - 3.1. EEG software viewing option, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
 - 3.2 Software option for viewing EEG, version: EASY (on request)
 - 3.3 Software option for viewing EEG in "off-line" mode, version: full, CFM, DSA (on request)

- 3.4 Software option for connecting to databases of external applications (optional)
- 3.5 EEG software option "Event identification" (as required)
- 3.6 EEG software option "Analysis of high-frequency oscillations" (on request)
- 3.7. Software option EEG "Brain Mapping", version: CLINIC, STUDIO (on request)
- 3.8. EEG software option "Data exchange with MRI" (on request)
- 3.9. Software EEG option "Interface HOLBERG SCORE" (optional)
- 3.10. Software EEG option "Audio playback" (optional)
- 3.11. Software EEG option "Averaging in "off-line" mode" (on demand)
- 3.12. Software option for viewing digital video-EEG, version: CLINIC, STUDIO (on request)
- 3.13. Software option for viewing digital video-EEG in Full HD format (on request)
- 3.14. Software option for viewing digital video-EEG in multichannel mode, version: CLINIC (on request)
- 3.15. Software option for advanced management of the archive server (optional)
- 3.16. Program statistics option (optional)
- 3.17. Software option for sleep program interface (optional)
- 3.18. Software option for background data transfer (optional)
- 3.19. Software option for exporting EEG data to the Viewer (optional)
- 3.20. Software option to enable HIS interface (optional)

4. Software options for sleep analysis:

- 4.1 Software option for analyzing imported curves in EDF format (optional)
- 4.2 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring SD LTM 32 EXPRESS (32 channels) (on request)
- 4.3 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring SD LTM 64 EXPRESS (64 channels) (optional)
- 4.4 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring 32 channels (on demand)
- 4.5 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring 64 channels (on request)
- 4.6 Software option for the analysis of curves recorded on the EEG Head Units of the SD PLUS series (on request)

5. Program control options.

- 5.1 Software option "Centralized control of video cameras" on a USB -key (on request)
- 5.2. Software option "Event management" on the USB -key (on demand)

6. USB key for registration and viewing programs (optional)

7. USB dongle for viewers for 1/5/10/25 users (optional)

8. USB dongle for sleep analysis programs (optional)

BQ interface module USB PLUS to EEG headunit for long-term SD monitoring LTM PLUS (64 channels)

10 Server type CPU

11. LCD display: 19"/22"/24"/27" (optional)

12. HBC EXPRESS cable for connecting an interface board or an interface module to EEG Head Units (optional)

13. EEG head unit SD PLUS FLEXI (38 channels)

15. BQ NET POE 2 interface module for SD PLUS EEG headunit connection, includes SUPP BQ NET POE 2 wall mount and BQ USB EXPRESS interface module
16. TERM EXPRESS quick-disconnect cable to the EEG Head Unit for long-term monitoring SD LTM EXPRESS in a configuration of up to 64 channels (on request)
17. BQ NET POWER POE power supply for BQ NET POE series interface modules (optional)
18. BQ PCI EXPR LTM FULL interface board based on PCI Express FULL Size to control SD PLUS EEG headunit
19. BQ PCI EXPR LTM HALF interface board based on PCI Express HALF Size to control SD PLUS EEG headunit
20. BQ USB EXPRESS Interface Module to SD PLUS EEG Head Unit
21. BOX TERMINAL EXPRESS interface module for SDPLUS EEG Head Units
22. USB adapter to serial port ADAPT USB SER (optional)
23. Isolated USB to Serial Adapter ADAPT USB SER (optional) ISO
24. CB MAG RAPID cable for connecting the Micromed system to Magstim magnetic stimulation devices (D-Sub 26 pin (optional)
25. Trigger cable CB RCRC4N for stimulator (optional)
26. Trigger cable CB RCRC4R for stimulator (optional)
27. Trigger cable CB RCRC1R for stimulator (optional)
28. Pair of CB RCRC5 PAIR trigger cables for stimulator (optional)
29. CB SER RCA EXT cable for synchronization with external stimulators (optional)
30. Fiber optic cable 10m/15m/20m/30m/40m/50m/100m (optional)
31. LED POWER power supply for LED photostimulator
32. 180 Extension cable LED EXT for LED photostimulator (on request)
33. CAMERA IP HD PRO video camera (optional)
34. CAMERA HD Full HD (1920 x 1080) camcorder with PTZ control and automatic IR switching (on demand)
35. CAMERA IP HD OUTDOOR video camera (optional)
36. CAMERA IP HD Full HD (1920 x 1080) camera with PTZ control and automatic IR switching for IP video (on demand)
37. Camcorder CAMERA LAPTOP with USB 2.0, including microphone (optional)
38. CAMERA IP COMPACT video camera (optional)
39. GC MM video card with dual monitor output (optional)
40. HDMI EXTENDER device for transmitting HDMI signal over Cat5e/6 cable, including receiver and transmitter (optional)
41. IP MIC microphone for IP video cameras (optional)
42. IP POWER power adapter for IP video cameras (optional)
43. VLAMP IR illuminator with auto-on (on demand)
44. Composite switcher VMIXER PRO with 2 inputs (optional)
45. Professional switcher VMIXER PRO with 5 HDMI inputs (optional)
46. ARM GM PLUS bracket for SD PLUS series EEG head units, complete
47. Vertical bracket ARM LCD DOUBLE for two LCD displays for fixing on the cart (optional)
48. Horizontal bracket ARM LCD HORIZ for two LCD displays for fixing on the cart (optional)
49. Vertical bracket ARM LCD SW for LCD displays for fixing on the trolley (optional)
50. HBC BT/RJ45 extension adapter for BOX TERMINAL EXPRESS interface module (module side)(as required)
51. Rack with base BASE STAND (optional)
52. BQ CLAMP (optional)
53. Swivel mount CAM JOIN for a video camera on a stand (optional)
54. CAM STAND for mounting a video camera or head unit on a trolley (optional)
55. FIX GM mount for trolley bracket (optional)
56. Fixing FIX GM TW for table bracket (optional)
57. JOIN GM connection for attaching the head unit holder to the bracket (optional)

58. HBC HT/RJ45 extension adapter for BOX TERMINAL EXPRESS interface module (PC side)(as required)
59. LED STAND PHOTIC wheel base for LED photostimulator (on request)
60. Bracket LED SUPP PHOTIC for fixing the LED photostimulator (optional)
61. Holder STAND CAM SW video camera assembly (optional)
62. Stand STAND POLE NEW on wheels with weighting (optional)
63. SUPP BQ NET POE2 wall mount for BQ NET POE 2 / BQ NET POE M interface modules (optional)
64. SUPP PLUS Holder for SD PLUS EEG Head Units (optional)
65. VCLAMP mount for video camera or IR light (optional)
66. VSTAND tripod stand for camcorder or IR illuminator (optional)
67. VSUPP Wall Mount for Video Camera (Optional)
68. VSUPP IP metal support for attaching an IP video camera to a swivel mount (optional)
69. VSUPP IP IR metal support for attaching the IP video camera to the rack (optional)
70. Wall mount VSUPP IP HD PLUS BLACK for video camera, black (optional)
71. VSUPP IP HD PLUS WHITE wall mount for video camera, white (optional)
72. ADAPT BNC RCA adapter (optional)
73. ADAPT DVI HDMI adapter (optional)
74. ANY WHERE 2 USB over LAN device with 2 USB ports (optional)
75. ANY WHERE 5 USB over LAN device with 5 USB ports (optional)
76. USB joystick CAM JOYSTICK for camera control (optional)
77. 24-key CAM PAD numeric keypad with USB connection for camera control (optional)
78. Numeric keypad CAM PAD 48 x 48 keys with USB connection for camera control (optional)
79. Connecting cable CB BBNB1 (BNC to BNC) (optional)
80. Connecting cable CB BNRC1 (BNC - RCA (Red) (as required)
81. Connecting cable CB BNRC5 (BNC - RCA (Black) (optional)
82. CB DVI cable for monitor (optional)
83. HDMI - HDMI MINI cable for video transmission, length 2m/5m/10m (optional)
84. HDMI to HDMI cable for high-definition video transmission, 2m/5m/15m long (optional)
85. Audio cable CB JACK JMINI with 6.3 mm plugs - MiniJack 3.5 mm (optional)
86. Audio cable CB JACK RCRC with 6.3 mm plugs - double RCA (optional)
87. Microphone cable CB MICROPHONE (optional)
88. CB MON cable for monitor (D-Sub 15 pin - D-Sub 15 pin) (optional)
89. Two-wire audio cable CB RCRC JMINI, RCA male-MiniJack 3.5 mm (optional)
90. Two-wire audio cable CB RCRC PAIR, RCA male-double RCA (as required)
91. CB VIDEO video cable as required)
92. CB IP HD PLUS 15 cable for controlling IP cameras (optional)
93. CB IP HD PLUS DC 5 cable for connecting IP HD PLUS cameras to each other (optional)
94. CB IP HD PLUS EXTENDER adapter for IP HD PLUS signal transmission over Cat5e/6 cable. pair (optional)
95. HDMI BOOSTER self-powered HDMI amplifier (optional)
96. 246 CARSTAR mobile trolley (on request)
97. 247 FIX TOR bracket for attaching the isolation transformer to the trolley (optional)
98. 248 Mouse pad small/regular/large (optional)
99. 249 Trolley CAR PLUS mobile, adjustable in height (as required)
100. FIX CART plates for fixing trolley legs (optional)
101. Power cable set KIT CB POWER for trolley (as required)
102. Set of labels KIT LABEL PLUS for the trolley (as required)
103. SUPP KEYB PLUS keyboard drawer for trolley (optional)
104. SUPP PC PLUS trolley mount (optional)
105. Small computer mount SUPP PC PLUS SFF on cart (optional)
106. Horizontal bar BAR POLE NEW for attaching the bracket to the trolley (optional)

107. Trolley CAR TORO with shelves, standard size (as required)
108. Trolley CAR TORO SMALL with shelves, small size (on request)
109. COVER TORO rear cover for standard CARTORO trolley (optional)
110. Back cover COVER TORO SMALL for small cart CAR TORO SMALL (as required)
111. Isolation transformer TOR 500, Pmax 500 VA, black (on request)
112. Isolation transformer TOR 500 W, Pmax 500 VA, white (on request)
113. Isolation transformer TR 50, 100 VA (on request)
114. NET ISOLAT network isolator for 10/100 Mbps networks, with RJ45 connectors (optional)
115. Network cable NET RJ45 with RJ45 connector, length 3m/5m/10m (optional)
116. 4-port network switch NET SWITCH 4 (optional)
117. 8-port network switch NET SWITCH 8 (optional)
118. Professional network switch NET SWITCH PRO (optional)
119. Standard Bluetooth USB interface (optional)
120. BS USB READER for Compact Flash card reader (optional)
121. PLIERS TP forceps for removing electrodes from the electrode box (optional)
122. USB extension cable USB EXT CB (optional)
123. USB EXTENDER device for transmitting a USB signal over a Cat5e / 6 cable, including a receiver and transmitter (optional)
124. BQ CONTACT contact closure control module (optional)
125. Patient button with 5m cable and RCA connector (optional)
126. BUTTON BOX to control patient buttons (optional)
127. Cables for BUTTON BOX (optional)
128. Illuminated patient button with 2 m cable and LEMO connector (optional)
129. Power supply for BUTTON BOX (optional)
130. Cable for connection unit for oximetric sensors. Length 3 m (as required)
131. Oximetry Sensor Connection Box for use with SD PLUS Series EEG Head Units (optional)
132. Oximetric sensor "EarClip", adult (on request)
133. FingerClip oximetry sensor, adult (optional)
134. Oximetric sensor "FingerClip", child (on request)
135. Fixing tape for SD OXY RA, adult size (as required)
136. Fixing tape SD OXY RC, child size (as required)
137. Flexible oximetry probe, adult size (optional)
138. Flexible oximetry sensor, child size (optional)
139. Large size reusable soft oximeter probe (digital readout 12.5 - 25.5 mm) (optional)
140. Reusable medium size soft oximeter probe (digital readout 10 - 19 mm) (optional)
141. Reusable soft oximetry sensor small size (digital display 7.5 - 12.5 mm) (optional)
142. KVM-extender COMP DVI RX compact, for receiving DVI + USB HDI (optional)
143. KVM Extender COMP DVI RX AUD compact, for receiving DVI + USB HDI + analog audio signal (optional)
144. KVM Extender COMP DVI TX compact, for DVI + USB HDI transmission (optional)
145. COMP DVI TX AUD compact KVM extender for DVI + USB HDI + analog audio transmission (optional)
146. DVI RX KVM module for DVI + USB HID reception (optional)
147. KVM Module DVI RX AUD Dual to receive DVI + USB HID + analog audio signal (as needed)
148. DVI TX KVM module for DVI + USB HID transmission (optional)
149. KVM Module DVI TX AUD Dual to Transmit DVI + USB HID + Analog Audio Signal (As Required)
150. KVM Module DVI TX AUD DCAT Dual to DVI + USB HID + Analog Audio Dual Channel Transmission (as needed)
151. DVI TX DCAT KVM Module for DVI + USB HID Dual Link Transmission (optional)

152. KIT KVM KVM module including transmitter and receiver for keyboard, mouse and monitor. (optional)
153. KVM module KIT KVM DUAL, including transmitter and receiver, for keyboard, mouse and two monitors (optional)
154. KVM Chassis KVM BODY 2 for 2 DVI modules (optional)
155. KVM Chassis KVM BODY 4 for 4 DVI modules (optional)
156. KVM BODY 4R Chassis for 4 DVI modules, including rack mount bracket (optional)
157. KVM Chassis KVM BODY 6 for 6 DVI modules (optional)
158. KVM BODY 6R Chassis for 6 DVI modules, including rack mount bracket (optional)
159. KVM IP Module, KVM IP RX, Receiver (optional)
160. KVM IP module, KVM IP TX, Transmitter (optional)
161. KVM matrix KVM SWITCH for 8/16/32/48/64/80 ports (on demand)

162. Electrode cap (optional):

- MRINEO 19 X
- MRINEO 31 X
- MRINEO 61 X

163. Cable for EL CFM electrodes (as required):

- CB CFM050
- CB CFM150

164. Electrodes with cable for LTM (as required):

- CB EL TPCO 19
- CB EL TPCO 31

165. Electrode (plate) EL CFM for 2-channel EEG (on demand)

166. Disk electrode with cable (as required):

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-EL TPCO SG6 200

-EL TPCO SG6 250

168. Ear electrode (clip) EL TPEAR with cable (on request)

169. EL TPLINK cable for combining two ear electrodes EL TPEAR (optional)

170. Connecting cable for LTM / LTM 256 (on request)

171. CB PUD150 cable for EL PUD electrodes (optional)

172. Disposable electrodes EL PUD (on demand)

173. Electrode monopolar needle (on demand):

-EMG MBOT2530

-EMG MBOT2540

-EMG MBOT3745

-EMG MBOT5050

-EMG MBOT7570

174. Snoring sensor (optional):

-EPM 1200

-EPM 1250

175. EPM 1312 breathing sensor kit, includes sensor (EPM 1387) and 4 tapes in different sizes (EPM 1340 + EPM 1341 + EPM 1342 + EPM 1343) (optional)

176. Pediatric breathing sensor EPM 1390 (on demand)

177. Pediatric air flow sensor EPM 1400 (on demand)

178. EPM 1439FM/ GM / HM Disposable Breath Sensor Kit, includes: EPM 1415 cable and 10 sensors for Adult 1426 EPMs / Medium Size 1427 EPMs / Infant Size 1428 EPMs) (as required)

179. Air flow sensor EPM 1450 (optional)

180. Set of air flow and pressure sensors EPM 1475 (as required)

181. Air pressure and snoring sensor EPM 14830 (on request)

182. Body position sensor EPM 1563 (optional)

183. EPM 1715 motion sensor kit (including piezoelectric transducer (EPM 1770) and pairs of ankle straps (EPMS 1730) (optional)

184. Chest force sensor kit, including interface cable and 4 chest straps (as required):

- EPM 9102 A for adults

- EPM 9102 C pediatric

185. Abdominal force sensor kit, including interface cable and 4 abdominal bands (as required):

- EPM 9101 A for adults

- EPM 9101 C pediatric

186. A set of abdominal and thoracic force sensors, including abdominal and thoracic interface cables, and 4 tapes (as required):

- EPM 914x A for adults

- EPM 914x C pediatric

- EPM 915x A for adults

- EPM 915x C pediatric

187. Tape for the breath sensor (as needed):

- EPMs 1340

- EPMs 1341

- EPMs 1342

- EPMs 1343

188. Piezocrystalline breath sensor EPMs 1387 (optional)

189. Cable for single-channel airflow sensor EPMs 1415 (optional)

190. Thermal breath sensor (on demand):

-EPMs 1426 adult size

-EPMs 1427 medium size

-EPMs 1428 pediatric

191. EPMS 1450-B interface for EPM 1450 airflow sensor (optional)

192. EPMS 1450- S transducer for EPM 1450 airflow sensor (optional)

193. Combined thermistor pressure and airflow sensor EPMs 1475

CO 2 cannula (as required):

- EPMs 14802, nasal and oral, male adult
- EPMs 14803, male for newborns
- EPMs 14804, male pediatric
- EPMs 14805, male adult

194. Piezoelectric pressure transducer for airflow and snoring EPMs 14830 (optional)

195. Elastic tape EPMs 1550 for body position sensor EPM 1520 (optional)

196. Converter EPMs 1563 for body position sensor EPM 1563 (optional)

197. ThermoCan cannula (as required):

- EPMs 15802

- EPMs 15807

198. Pair of EPMs 1730 ankle straps (optional)

199. Piezoelectric transducer EPMs 1770 for motion sensor EPM 1715 (optional)

200. A pair of induction tapes (as needed):

- 9002-L120 EPMs

- 9002-L133 EPMs

- 9002-L40 EPMs

- 9002-L60 EPMs

- 9002-L90 EPMs

201. Adjustable induction tape EPMs 9010 (as required)

202. Induction interface cable (as required):

- EPMs 9101 for the abdomen

- EPMs 9102 for the thoracic region

203. Disposable ECG electrode ECG EL (optional)

204. Cable for ECG electrode (optional):

-ECG EL TPCLB black

-ECG EL TPCLR red

VII. ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL STUDIES MICROMED , PERFORMED BY: SD PLUS FLEXI HIGH RATE , IN THE COMPOSITION

System software PLUS EVOLUTION - BRAIN QUICK including the following software options:

Software option for registration and viewing:

2.1 Software option for recording and viewing EEG, 2-channel CFM version (on request)

2.2 Software option for recording and viewing EEG, 8-channel CFM version (on request)

2.3 Software option for recording and viewing EEG, version: CLINIC, STUDIO, CLINIC MULTICH (on request)

2.4 Software option for recording and viewing digital video-EEG, version CLINI C (on request)

2.5. Software option for recording and viewing digital video-EEG in FullHD format (on request)

2.6 Software option for recording and viewing digital video-EEG over the network, version CLINI C (on request)

2.7 Software option for real-time logging and viewing over the network (as required)

2.8 Software option for EEG analysis in "on-line" and "off-line" modes, version: full / CFM / DSA. (on demand)

2.9 Software option for transmitting EEG data via TCP/IP (optional)

2.10 Software option " ERP -manager" (optional)

2.11 Software option "Management of outpatient registrar" (as required)

3. Software options for viewing:

- 3.1. EEG software viewing option, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
- 3.2 Software option for viewing EEG, version: EASY (on request)
- 3.3 Software option for viewing EEG in "off-line" mode, version: full, CFM, DSA (on request)
- 3.4 Software option for connecting to databases of external applications (optional)
- 3.5 EEG software option "Event identification" (as required)
- 3.6 EEG software option "Analysis of high-frequency oscillations" (on request)
- 3.7. Software option EEG "Brain Mapping", version: CLINIC, STUDIO (on request)
- 3.8. EEG software option "Data exchange with MRI" (on request)
- 3.9. Software EEG option "Interface HOLBERG SCORE" (optional)
- 3.10. Software EEG option "Audio playback" (optional)
- 3.11. Software EEG option "Averaging in "off-line" mode" (on demand)
- 3.12. Software option for viewing digital video-EEG, version: CLINIC, STUDIO (on request)
- 3.13. Software option for viewing digital video-EEG in Full HD format (on request)
- 3.14. Software option for viewing digital video-EEG in multichannel mode, version: CLINIC (on request)
- 3.15. Software option for advanced management of the archive server (optional)
- 3.16. Program statistics option (optional)
- 3.17. Software option for sleep program interface (optional)
- 3.18. Software option for background data transfer (optional)
- 3.19. Software option for exporting EEG data to the Viewer (optional)
- 3.20. Software option to enable HIS interface (optional)

4. Software options for sleep analysis:

- 4.1 Software option for analyzing imported curves in EDF format (optional)
- 4.2 Software option for the analysis of curves recorded on the EEG Head Units of the SD PLUS series (on request)

5. Program control options.

- 5.1 Software option "Centralized control of video cameras" on a USB -key (on request)
 - b. Software option "Event management" on a USB key (optional)

6. USB key for registration and viewing programs (optional)

7. USB dongle for viewers for 1/5/10/25 users (optional)

8. USB dongle for sleep analysis programs (optional)

BQ interface module USB PLUS to EEG head unit for long-term monitoring SDLTPLUS (64 channels)

10 HBC EXPRESS cable for connecting an interface board or interface module to EEG Head Units (optional)

11. LCD display: 19"/22"/24"/27" (optional)

12. TERM EXPRESS quick-disconnect cable to the EEG Head Unit for long-term monitoring SD LTM EXPRESS in a configuration of up to 64 channels (on request)

13. EEG head unit SD PLUS FLEXI HIGH RATE (29 channels)

14. HBC BT/RJ45 extension adapter for BOX TERMINAL EXPRESS interface module (module side)(as required)

15. Photostimulator LED with bracket

16. BQ NET POE 2 interface module for SD PLUS EEG headunit connection, includes SUPP BQ NET POE 2 wall mount and BQ USB EXPRESS interface module
17. - HBC HT/RJ45 extension adapter for BOX TERMINAL EXPRESS interface module (PC side) (optional)
18. BQ NET POWER POE power supply for BQ NET POE series interface modules (optional)
19. BQ PCI EXPR LTM FULL interface board based on PCI Express FULL Size to control SD PLUS EEG headunit
20. BQ PCI EXPR LTM HALF interface board based on PCI Express HALF Size to control SD PLUS EEG headunit
21. BQ USB EXPRESS Interface Module to SD PLUS EEG Head Unit
22. BOX TERMINAL EXPRESS interface module for SDPLUS EEG Head Units
23. USB adapter to serial port ADAPT USB SER (optional)
24. Isolated USB to Serial Adapter ADAPT USB SER (optional) ISO
25. CB MAG RAPID cable for connecting the Micromed system to Magstim magnetic stimulation devices (D-Sub 26 pin) (optional)
26. Trigger cable CB RCRC4N for stimulator (optional)
27. Trigger cable CB RCRC4R for stimulator (optional)
28. Trigger cable CB RCRC1R for stimulator (optional)
29. Pair of CB RCRC5 PAIR trigger cables for stimulator (optional)
30. CB SER RCA EXT cable for synchronization with external stimulators
31. Fiber optic cable 10m/15m/20m/30m/40m/50m/100m (optional)
32. LED POWER power supply for LED photostimulator
33. 180 Extension cable LED EXT for LED photostimulator (on request)
34. CAMERA IP HD PRO video camera (optional)
35. CAMERA HD Full HD (1920 x 1080) camcorder with PTZ control and automatic IR switching (on demand)
36. CAMERA IP HD OUTDOOR video camera (optional)
37. CAMERA IP HD PLUS video camera (optional)
38. Camcorder CAMERA LAPTOP with USB 2.0, including microphone (optional)
39. CAMERA IP COMPACT video camera (optional)
40. GC MM video card with dual monitor output (optional)
41. HDMI EXTENDER device for transmitting HDMI signal over Cat5e/6 cable, including receiver and transmitter (optional)
42. IP MIC microphone for IP video cameras (optional)
43. IP POWER power adapter for IP video cameras (optional)
44. VLAMP IR illuminator with auto-on (on demand)
45. Composite switcher VMIXER PRO with 2 inputs (optional)
46. Professional switcher VMIXER PRO with 5 HDMI inputs (optional)
47. ARM GM PLUS bracket for SD PLUS series EEG head units. complete
48. Vertical bracket ARM LCD DOUBLE for two LCD displays for fixing on the cart (optional)
49. Horizontal bracket ARM LCD HORIZ for two LCD displays for fixing on the cart (optional)
50. Vertical bracket ARM LCD SW for LCD displays for fixing on the trolley (optional)
51. Rack with base BASE STAND (optional)
52. BQ CLAMP (optional)
53. Swivel mount CAM JOIN for a video camera on a stand (optional)
54. CAM STAND for mounting a video camera or head unit on a trolley (optional)
55. FIX GM mount for trolley bracket (optional)
56. Fixing FIX GM TW for table bracket (optional)
57. JOIN GM connection for attaching the head unit holder to the bracket (optional)
58. LED STAND PHOTIC wheel base for LED photostimulator (on demand)
59. Bracket LED SUPP PHOTIC for fixing the LED photostimulator (optional)
60. Holder STAND CAM SW video camera assembly (optional)

61. Stand STAND POLE NEW on wheels with weighting (optional)
62. SUPP BQ NET POE2 wall mount for BQ NET POE 2 interface modules BQ NET POE 2 / BQ NET POE M (optional) (optional)
63. SUPP PLUS Holder for SD PLUS EEG Head Units (optional)
64. VCLAMP mount for video camera or IR light (optional)
65. VSTAND tripod stand for camcorder or IR illuminator (optional)
66. VSUPP Wall Mount for Video Camera (Optional)
67. VSUPP IP metal support for attaching an IP video camera to a swivel mount (optional)
68. VSUPP IP IR metal support for attaching the IP video camera to the rack (optional)
69. Wall mount VSUPP IP HD PLUS BLACK for video camera, black (optional)
70. VSUPP IP HD PLUS WHITE wall mount for video camera, white (optional)
71. ADAPT BNC RCA adapter (optional)
72. ADAPT DVI HDMI adapter (optional)
73. ANY WHERE 2 USB over LAN device with 2 USB ports (optional)
74. ANY WHERE 5 USB over LAN device with 5 USB ports (optional)
75. USB joystick CAM JOYSTICK for camera control (optional)
76. 24-key CAM PAD numeric keypad with USB connection for camera control (optional)
77. Numeric keypad CAM PAD 48 x 48 keys with USB connection for camera control (optional)
78. Connecting cable CB BBNB1 (BNC to BNC) (optional)
79. Connecting cable CB BNRC1 (BNC - RCA (Red) (as required)
80. Connecting cable CB BNRC5 (BNC - RCA (Black) (optional)
81. CB DVI cable for monitor (optional)
82. HDMI - HDMI MINI cable for video transmission, length 2m/5m/10m (optional)
83. HDMI to HDMI cable for high-definition video transmission, 2m/5m/15m long (optional)
84. Audio cable CB JACK JMINI with 6.3 mm plugs - MiniJack 3.5 mm (optional)
85. Audio cable CB JACK RCRC with 6.3 mm plugs - double RCA (optional)
86. Microphone cable CB MICROPHONE (optional)
87. CB MON cable for monitor (D-Sub 15 pin - D-Sub 15 pin) (optional)
88. Two-wire audio cable CB RCRC JMINI, RCA male-MiniJack 3.5 mm (optional)
89. Two-wire audio cable CB RCRC PAIR, RCA male-double RCA (as required)
90. CB VIDEO video cable (optional)
91. CB IP HD PLUS 15 cable for controlling IP cameras (optional)
92. CB IP HD PLUS DC 5 cable for connecting IP HD PLUS cameras to each other (optional)
93. CB IP HD PLUS EXTENDER adapter for IP HD PLUS signal transmission over Cat5e/6 cable, pair (optional)
94. HDMI BOOSTER self-powered HDMI amplifier (optional)
95. Trolley CARSTAR mobile (on demand)
96. Bracket FIX TOR for attaching the isolation transformer to the trolley (optional)
97. Mouse pad small/regular/large (as needed)
98. Trolley CAR PLUS mobile, adjustable in height (as required)
99. FIX CART plates for fixing trolley legs (optional)
100. Power cable set KIT CB POWER for trolley (as required)
101. Set of labels KIT LABEL PLUS for the trolley (as required)
102. SUPP KEYB PLUS keyboard drawer for trolley (optional)
103. SUPP PC PLUS trolley mount (optional)
104. Small computer mount SUPP PC PLUS SFF on cart (optional)
105. Horizontal bar BAR POLE NEW for attaching the bracket to the trolley (optional)
106. Trolley CAR TORO with shelves, standard size (as required)
107. Trolley CAR TORO SMALL with shelves, small size (on demand)
108. COVER TORO rear cover for standard CARTORO trolley (optional)
109. Back cover COVER TORO SMALL for small cart CAR TORO SMALL (as required)
110. Isolation transformer TOR 500, Pmax 500 VA, black (on request)

111. Isolation transformer TOR 500 W, Pmax 500 VA, white (on request)
112. Isolation transformer TR 50, 100 VA (on request)
113. NET ISOLAT network isolator for 10/100 Mbps networks, with RJ45 connectors (optional)
114. Network cable NET RJ45 with RJ45 connector, length 3m/5m/10m (optional)
115. 4-port network switch NET SWITCH 4 (optional)
116. 8-port network switch NET SWITCH 8 (optional)
117. Professional network switch NET SWITCH PRO (optional)
118. Standard Bluetooth USB interface (optional)
119. BS USB READER for Compact Flash card reader (optional)
120. PLIERS TP forceps for removing electrodes from the electrode box (optional)
121. USB extension cable USB EXT CB (optional)
122. USB EXTENDER device for transmitting a USB signal over a Cat5e / 6 cable, including a receiver and transmitter (optional)
123. BQ CONTACT contact closure control module (optional)
124. Patient button with 5m cable and RCA connector (optional)
125. BUTTON BOX to control patient buttons (optional)
126. Cables for BUTTON BOX (optional)
127. Illuminated patient button with 2 m cable and LEMO connector (optional)
128. Power supply for BUTTON BOX (optional)
129. Cable for connection unit for oximetric sensors. Length 3 m (as required)
130. Oximetry Sensor Connection Box for use with SD PLUS Series EEG Head Units (optional)
131. Oximetric sensor "EarClip", adult (on request)
132. FingerClip oximetry sensor, adult (optional)
133. Oximetric sensor "FingerClip", child (on request)
134. Fixing tape for SD OXY RA, adult size (as required)
135. Fixing tape SD OXY RC, child size (as required)
136. Flexible oximetry probe, adult size (optional)
137. Flexible oximetry sensor, child size (optional)
138. Large size reusable soft oximeter probe (digital readout 12.5 - 25.5 mm) (optional)
139. Reusable medium size soft oximeter probe (digital readout 10 - 19 mm) (optional)
140. Reusable soft oximetry sensor small size (digital display 7.5 - 12.5 mm) (optional)
141. KVM-extender COMP DVI RX compact, for receiving DVI + USB HDI (optional)
142. KVM Extender COMP DVI RX AUD compact, for receiving DVI + USB HDI + analog audio signal (optional)
143. KVM Extender COMP DVI TX compact, for DVI + USB HDI transmission (optional)
144. COMP DVI TX AUD compact KVM extender for DVI + USB HDI + analog audio transmission (optional)
145. DVI RX KVM module for DVI + USB HID reception (optional)
146. KVM Module DVI RX AUD Dual to receive DVI + USB HID + analog audio signal (as required)
147. DVI TX KVM module for DVI + USB HID transmission (optional)
148. KVM Module DVI TX AUD Dual to Transmit DVI + USB HID + Analog Audio Signal (As Required)
149. KVM Module DVI TX AUD DCAT Dual to DVI + USB HID + Analog Audio Dual Channel Transmission (as required)
150. DVI TX DCAT KVM Module for DVI + USB HID Transmission, with Dual Channel Transmission (optional)
151. KIT KVM KVM module including transmitter and receiver for keyboard, mouse and monitor. (optional)
152. KVM module KIT KVM DUAL, including transmitter and receiver, for keyboard, mouse and two monitors (optional)
153. KVM Chassis KVM BODY 2 for 2 DVI modules (optional)

- 154.KVM Chassis KVM BODY 4 for 4 DVI modules (optional)
- 155.KVM BODY 4R Chassis for 4 DVI modules, including rack mount bracket (optional)
- 156.KVM Chassis KVM BODY 6 for 6 DVI modules (optional)
- 157.KVM BODY 6R Chassis for 6 DVI modules, including rack mount bracket (optional)
- 158.KVM IP Module, KVM IP RX, Receiver (optional)
- 159.KVM IP module, KVM IP TX, Transmitter (optional)
- 160.KVM matrix KVM SWITCH for 8/16/32/48/64/80 ports (on demand)
- 161. Electrode cap (optional):
 - MRINEO 19 X
 - MRINEO 31 X
 - MRINEO 61 X
- 161.Cable for EL CFM electrodes (as required):
 - CB CFM050
 - CB CFM150
- 162. Electrodes with cable for LTM (as required):
 - CB EL TPCO 19
 - CB EL TPCO 31
- 163. Electrode (plate) EL CFM for 2-channel EEG (on demand)
- 164. Disk electrode with cable (as required):
 - EL TPCO
 - EL TPCO 6
 - EL TPCO G06
 - EL TPCO G10
 - EL TPCO S 100
 - EL TPCO S 150
 - EL TPCO S 200
 - EL TPCO S 250
 - EL TPCO SG 100
 - EL TPCO SG 150
 - EL TPCO SG 200
 - EL TPCO SG 250
 - EL TPCO SLT 100
 - EL TPCO SLT 150
 - EL TPCO SLT 200
 - EL TPCO SLT 250
 - EL TPCO SLTG 100
 - EL TPCO SLTG 150
 - EL TPCO SLTG 200
 - EL _ TPCO SLTG 250
- 165. Pediatric disk electrode with cable (as required):
 - EL TPCO S6 100
 - EL TPCO S6 150
 - EL TPCO S6 200
 - EL TPCO S6 250
 - EL TPCO SG6 100
 - EL TPCO SG6 150
 - EL TPCO SG6 200
 - EL TPCO SG6 250
- 166. Ear electrode (clip) EL TPEAR with cable (on request)
- 167. EL TPLINK cable for combining two ear electrodes EL TPEAR (optional)
- 168.Connecting cable for LTM / LTM 256 (on request)
- 169.CB PUD150 cable for EL PUD electrodes (optional)

170. Disposable electrodes EL PUD (on demand)

171. Electrode monopolar needle (on demand):

- EMG MBOT2530

- EMG MBOT2540

- EMG MBOT3745

- EMG MBOT5050

- EMG MBOT7570

172. Snoring sensor (on demand):

- EPM 1200

- EPM 1250

173. EPM 1312 breathing sensor kit, includes sensor (EPM 1387) and 4 tapes in different sizes (EPM 1340 + EPM 1341 + EPM 1342 + EPM 1343) (optional)

174. Pediatric breathing sensor EPM 1390 (on demand)

175. Pediatric air flow sensor EPM 1400 (on request)

176. Set of EPM 1439FM/ GM / HM disposable breath sensors , includes: EPM 1415 cable and 10 sensors adult 1426 EPMs / medium size 1427 EPMs / infant size 1428 EPMs) (as required)

177. Air flow sensor EPM 1450 (optional)

178. Set of air flow and pressure sensors EPM 1475 (as required)

179. Air pressure sensor and snoring EPM 14830 (optional)

180. Body position sensor EPM 1563 (optional)

181. EPM 1715 motion sensor kit (including piezoelectric transducer (EPM 1770) and pairs of ankle straps (EPMS 1730) (optional)

182. Chest force sensor kit, including interface cable and 4 chest straps (as required):

- EPM 9102 A for adults

- EPM 9102 C pediatric

183. Set of abdominal force sensor, including interface cable and 4 abdominal bands (as required):

- EPM 9101 A for adults

- EPM 9101 C pediatric

184. A set of abdominal and thoracic force sensors, including abdominal and thoracic interface cables, and 4 tapes (as required):

- EPM 914x A for adults

- EPM 914x C pediatric

- EPM 915x A for adults

- EPM 915x C pediatric

185. Tape for the breath sensor (as needed):

- EPMs 1340

- EPMs 1341

- EPMs 1342

- EPMs 1343

186. Piezocrystalline breath sensor EPMs 1387 (optional)

187.. Cable for single-channel airflow sensor EPMs 1415 (on request)

188. Thermal breath sensor (as required):

- EPMs 1426 adult size

- EPMs 1427 medium size

- EPMs 1428 pediatric

189. EPMs 1450-B interface for EPM 1450 air flow sensor (optional)

S transducer for EPM 1450 airflow sensor (optional)

191. Combined thermistor pressure and air flow sensor EPMs 1475

192. CO 2 cannula (as required):

- EPMs 14802, nasal and oral, male adult

- EPMs 14803, male for newborns

- EPMs 14804, male pediatric

- EPMs 14805, male adult

- 193. Piezoelectric pressure transducer for air flow and snoring EPMs 14830 (on request)
- 194. Elastic tape EPMs 1550 for body position sensor EPM 1520 (optional)
- 195. EPMs 1563 transducer for EPM 1563 body position sensor (optional)
- 196. ThermoCan cannula (optional):
 - EPMs 15802
 - EPMs 15807
- 197. Pair of EPMs 1730 ankle straps (optional)
- 198. Piezoelectric transducer EPMs 1770 for motion sensor EPM 1715 (optional)
- 199. A pair of induction tapes (as needed):
 - 9002-L120 EPMs
 - 9002-L133 EPMs
 - 9002-L40 EPMs
 - 9002-L60 EPMs
 - 9002-L90 EPMs
- 200. Adjustable induction band EPMs 9010 (as required)
- 201. Induction interface cable (as required):
 - EPMs 9101 for the abdomen
 - EPMs 9102 for the thoracic region
- 202. Disposable ECG electrode ECG EL (on demand)
- 203.. ECG electrode cable (optional):
 - ECG EL TPCLB black
 - ECG EL TPCLR red

VIII. ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL STUDIES MICROMED , PERFORMED BY: SD PLUS RESEARCH, INCLUDED

System software PLUS EVOLUTION - BRAIN QUICK including the following software options:

Software option for registration and viewing:

- 2.1 Software option for recording and viewing EEG, 2-channel CFM version (on request)
- 2.2 Software option for recording and viewing EEG, 8-channel CFM version (on request)
- 2.3 Software option for recording and viewing EEG, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
- 2.4 Software option for recording and viewing digital video-EEG, version CLINI C (on request)
- 2.5. Software option for recording and viewing digital video-EEG in Full HD format (on request)
- 2.6 Software option for recording and viewing digital video-EEG over the network, version CLINI C (on request)
- 2.7 Software option for real-time logging and viewing over the network (as required)
- 2.8 Software option for EEG analysis in "on-line" and "off-line" modes, version: full / CFM / DSA. (on demand)
- 2.9 Software option for transmitting EEG data via TCP/IP (optional)
- 2.10 Software option " ERP -manager" (optional)
- 2.11 Software option "Management of outpatient registrar" (as required)
- 3. Software options for viewing:
 - 3.1. EEG software viewing option, version: CLINIC, STUDIO, CLINIC MULTICH (on request)

- 3.3 Software option for viewing EEG in "off-line" mode, version: full, CFM, DSA (on request)
 - 3.4 Software option for connecting to databases of external applications (optional)
 - 3.5 EEG software option "Event identification" (as required)
 - 3.6 EEG software option "Analysis of high-frequency oscillations" (on request)
 - 3.7. Software option EEG "Brain Mapping", version: CLINIC, STUDIO (on request)
 - 3.8. EEG software option "Data exchange with MRI" (on request)
 - 3.9. Software EEG option "Interface IOLBERG SCORE" (optional)
 - 3.10. Software EEG option "Audio playback" (optional)
 - 3.11. Software EEG option "Averaging in "off-line" mode" (on demand)
 - 3.12. Software option for viewing digital video-EEG, version: CLINIC, STUDIO (on request)
 - 3.13. Software option for viewing digital video-EEG in Full HD format (on request)
 - 3.14. Software option for viewing digital video-EEG in multichannel mode, version: CLINIC (on request)
 - 3.15. Software option for advanced management of the archive server (optional)
 - 3.16. Program statistics option (optional)
 - 3.17. Software option for sleep program interface (optional)
 - 3.18. Software option for background data transfer (optional)
 - 3.19. Software option for exporting EEG data to the Viewer (optional)
 - 3.20. Software option to enable HIS interface (optional)
4. Software options for sleep analysis:
- 4.1 Software option for analyzing imported curves in EDF format (optional)
 - 4.2 Software option for the analysis of curves recorded on the EEG Head Units of the SD PLUS series (on request)
5. Program control options.
- 5.1 Software option "Centralized control of video cameras" on a USB -key (on request)
 - c. Software option "Event management" on a USB key (optional)
6. USB key for registration and viewing programs (optional)
7. USB dongle for viewers for 1/5/10/25 users (optional)
8. USB dongle for sleep analysis programs (optional)
- BQUSBPLUS interface module to EEG head unit for long-term monitoring SDLTPLUS (64 channels)

10 HBC EXPRESS cable for connecting an interface board or interface module to EEG Head Units (optional)

11. LCD display: 19"/22"/24"/27" (optional)

12. EEG head unit SD PLUS RESEARCH (48 channels)

13. Photostimulator LED with bracket

14. BQ NET POE 2 interface module for SD PLUS EEG headunit connection, includes SUPP BQ NET POE 2 wall mount and BQ USB EXPRESS interface module

15. BQ NET POWER POE power supply for BQ NET POE series interface modules (optional)

16. BQ PCI EXPR LTM FULL interface board based on PCI Express FULL Size to control SD PLUS EEG headunit

17. BQ PCI EXPR LTM HALF interface board based on PCI Express HALF Size to control SD

18. BQ USB EXPRESS Interface Module to SD PLUS EEG Head Unit
19. BOX TERMINAL EXPRESS interface module for SD PLUS EEG Head Units
20. TERM EXPRESS quick-disconnect cable to the EEG head unit for long-term monitoring SD LTM EXPRESS in a configuration of up to 64 channels (on request)
21. HBC BT/RJ45 extension adapter for BOX TERMINAL EXPRESS interface module (module side) (optional)
22. HBC HT/RJ45 extension adapter for BOX TERMINAL EXPRESS interface module (PC side) (optional)
23. USB adapter to serial port ADAPT USB SER (optional)
24. Isolated USB to Serial Adapter ADAPT USB SER (optional) ISO
25. CB MAG RAPID cable for connecting the Micromed system to Magstim magnetic stimulation devices (D-Sub 26 pin (optional)
26. Trigger cable CB RCRC4N for stimulator (optional)
27. Trigger cable CB RCRC4R for stimulator (optional)
28. Trigger cable CB RCRC1R for stimulator (optional)
29. Pair of CB RCRC5 PAIR trigger cables for stimulator (optional)
30. CB SER RCA EXT cable for synchronization with external stimulators
31. Fiber optic cable 10m/15m/20m/30m/40m/50m/100m (optional)
32. LED POWER power supply for LED photostimulator
33. 180 Extension cable LED EXT for LED photostimulator (on request)
34. CAMERA IP HD PRO video camera (optional)
35. CAMERA HD Full HD (1920 x 1080) camcorder with PTZ control and automatic IR switching (on demand)
36. CAMERA IP HD OUTDOOR video camera (optional)
37. CAMERA IP HD PLUS video camera (optional)
38. Camcorder CAMERA LAPTOP with USB 2.0, including microphone (optional)
39. CAMERA IP COMPACT video camera (optional)
40. GC MM video card with dual monitor output (optional)
41. HDMI EXTENDER device for transmitting HDMI signal over Cat5e/6 cable, including receiver and transmitter (optional)
42. IP MIC microphone for IP video cameras (optional)
43. IP POWER power adapter for IP video cameras (optional)
44. VLAMP IR illuminator with auto-on (on demand)
45. Composite switcher VMIXER PRO with 2 inputs (optional)
46. Professional switcher VMIXER PRO with 5 HDMI inputs (optional)
47. ARM GM PLUS bracket for SD PLUS series EEG head units, complete
48. Vertical bracket ARM LCD DOUBLE for two LCD displays for fixing on the cart (optional)
49. Horizontal bracket ARM LCD HORIZ for two LCD displays for fixing on the cart (optional)
50. Vertical bracket ARM LCD SW for LCD displays for fixing on the trolley (optional)
51. Rack with base BASE STAND (optional)
52. BQ CLAMP (optional)
53. Swivel mount CAM JOIN for a video camera on a stand (optional)
54. CAM STAND for mounting a video camera or head unit on a trolley (optional)
55. FIX GM mount for trolley bracket (optional)
56. Fixing FIX GM TW for table bracket (optional)
57. JOIN GM connection for attaching the head unit holder to the bracket (optional)
58. LED STAND PHOTIC wheel base for LED photostimulator (on demand)
59. Bracket LED SUPP PHOTIC for fixing the LED photostimulator (optional)
60. Holder STAND CAM SW video camera assembly (optional)
61. Stand STAND POLE NEW on wheels with weighting (optional)
62. SUPP BQ NET POE2 wall mount for BQ NET POE 2 / BQ NET POE M interface modules (optional)

(optional)

63. SUPP PLUS Holder for SD PLUS EEG Head Units (optional)
64. VCLAMP mount for video camera or IR light (optional)
65. VSTAND tripod stand for camcorder or IR illuminator (optional)
66. VSUPP Wall Mount for Video Camera (Optional)
67. VSUPP IP metal support for attaching an IP video camera to a swivel mount (optional)
68. VSUPP IP IR metal support for attaching the IP video camera to the rack (optional)
69. Wall mount VSUPP IP HD PLUS BLACK for video camera, black (optional)
70. VSUPP IP HD PLUS WHITE wall mount for video camera, white (optional)
71. ADAPT BNC RCA adapter (optional)
72. ADAPT DVI HDMI adapter (optional)
73. ANY WHERE 2 USB over LAN device with 2 USB ports (optional)
74. ANY WHERE 5 USB over LAN device with 5 USB ports (optional)
75. USB joystick CAM JOYSTICK for camera control (optional)
76. 24-key CAM PAD numeric keypad with USB connection for camera control (optional)
77. Numeric keypad CAM PAD 48 x 48 keys with USB connection for camera control (optional)
78. Connecting cable CB BBN1 (BNC to BNC) (optional)
79. Connecting cable CB BNRC1 (BNC - RCA (Red) (as required)
80. Connecting cable CB BNRC5 (BNC - RCA (Black) (optional)
81. CB DVI cable for monitor (optional)
82. HDMI - HDMI MINI cable for video transmission. length 2m/5m/10m (optional)
83. HDMI to HDMI cable for high-definition video transmission, 2m/5m/15m long (optional)
84. Audio cable CB JACK JMINI with 6.3 mm plugs - MiniJack 3.5 mm (optional)
85. Audio cable CB JACK RCRC with 6.3 mm plugs - double RCA (optional)
86. Microphone cable CB MICROPHONE (optional)
87. CB MON cable for monitor (D-Sub 15 pin - D-Sub 15 pin) (optional)
88. Two-wire audio cable CB RCRC JMINI, RCA male-MiniJack 3.5 mm (optional)
89. Two-wire audio cable CB RCRC PAIR, RCA male-double RCA (as required)
90. CB VIDEO video cable as required)
91. CB IP HD PLUS 15 cable for controlling IP cameras (optional)
92. CB IP HD PLUS DC 5 cable for connecting IP HD PLUS cameras to each other (optional)
93. CB IP HD PLUS EXTENDER adapter for IP HD PLUS signal transmission over Cat5e/6 cable, pair (optional)
94. HDMI BOOSTER self-powered HDMI amplifier (optional)
95. Trolley CARSTAR mobile (on demand)
96. Bracket FIX TOR for attaching the isolation transformer to the trolley (optional)
97. Mouse pad small/regular/large (as needed)
98. Trolley CAR PLUS mobile, adjustable in height (as required)
99. FIX CART plates for fixing trolley legs (optional)
100. Power cable set KIT CB POWER for trolley (as required)
101. Set of labels KIT LABEL PLUS for the trolley (as required)
102. SUPP KEYB PLUS keyboard drawer for trolley (optional)
103. SUPP PC PLUS trolley mount (optional)
104. Small computer mount SUPP PC PLUS SFF on cart (optional)
105. Horizontal bar BAR POLE NEW for attaching the bracket to the trolley (optional)
106. Trolley CAR TORO with shelves, standard size (as required)
107. Trolley CAR TORO SMALL with shelves, small size (on demand)
108. COVER TORO rear cover for standard CARTORO trolley (optional)
109. Back cover COVER TORO SMALL for small cart CAR TORO SMALL (as required)
110. Isolation transformer TOR 500, Pmax 500 VA, black (on request)
111. Isolation transformer TOR 500 W, Pmax 500 VA, white (on request)
112. Isolation transformer TR 50, 100 VA (on request)

113. NET ISOLAT network isolator for 10/100 Mbps networks, with RJ45 connectors (optional)
114. Network cable NET RJ45 with RJ45 connector, length 3m/5m/10m (optional)
115. 4-port network switch NET SWITCH 4 (optional)
116. 8-port network switch NET SWITCH 8 (optional)
117. Professional network switch NET SWITCH PRO (optional)
118. Standard Bluetooth USB interface (optional)
119. BS USB READER for Compact Flash card reader (optional)
120. PLIERS TP forceps for removing electrodes from the electrode box (optional)
121. USB extension cable USB EXT CB (optional)
122. USB EXTENDER device for transmitting a USB signal over a Cat5e / 6 cable, including a receiver and transmitter (optional)
123. BQ CONTACT contact closure control module (optional)
124. Patient button with 5m cable and RCA connector (optional)
125. BUTTON BOX to control patient buttons (optional)
126. Cables for BUTTON BOX (optional)
127. Illuminated patient button with 2 m cable and LEMO connector (optional)
128. Power supply for BUTTON BOX (optional)
129. Cable for connection unit for oximetric sensors. Length 3 m (as required)
130. Oximetry Sensor Connection Box for use with SD PLUS Series EEG Head Units (optional)
131. Oximetric sensor "EarClip", adult (on request)
132. FingerClip oximetry sensor, adult (optional)
133. Oximetric sensor "FingerClip", child (on request)
134. Fixing tape for SD OXY RA, adult size (as required)
135. Fixing tape SD OXY RC, child size (as required)
136. Flexible oximetry probe, adult size (optional)
137. Flexible oximetry sensor, child size (optional)
138. Large size reusable soft oximeter probe (digital readout 12.5 - 25.5 mm) (optional)
139. Reusable medium size soft oximeter probe (digital readout 10 - 19 mm) (optional)
140. Reusable soft oximetry sensor small size (digital display 7.5 - 12.5 mm) (optional)
141. KVM-extender COMP DVI RX compact, for receiving DVI + USB HDI (optional)
142. KVM Extender COMP DVI RX AUD compact, for receiving DVI + USB HID + analog audio signal (optional)
143. KVM Extender COMP DVI TX compact, for DVI + USB HID transmission (optional)
144. COMP DVI TX AUD compact KVM extender for DVI + USB HID + analog audio transmission (optional)
145. DVI RX KVM module for DVI + USB HID reception (optional)
146. KVM Module DVI RX AUD Dual to receive DVI + USB HID + analog audio signal (as required)
147. DVI TX KVM module for DVI + USB HID transmission (optional)
148. KVM Module DVI TX AUD Dual to Transmit DVI + USB HID + Analog Audio Signal (As Required)
149. KVM Module DVI TX AUD DCAT Dual to DVI + USB HID + Analog Audio Dual Channel Transmission (as required)
150. DVI TX DCAT KVM Module for DVI + USB HID Dual Link Transmission
151. KIT KVM KVM module including transmitter and receiver for keyboard, mouse and monitor. (optional)
152. KVM module KIT KVM DUAL, including transmitter and receiver, for keyboard, mouse and two monitors (optional)
153. KVM Chassis KVM BODY 2 for 2 DVI modules (optional)
154. KVM Chassis KVM BODY 4 for 4 DVI modules (optional)

155. KVM BODY 4R Chassis for 4 DVI modules, including rack mount bracket (optional)
156. KVM Chassis KVM BODY 6 for 6 DVI modules (optional)
157. KVM BODY 6R Chassis for 6 DVI modules, including rack mount bracket (optional)
158. KVM IP Module, KVM IP RX, Receiver (optional)
159. KVM IP module, KVM IP TX, Transmitter (optional)
160. KVM matrix KVM SWITCH for 8/16/32/48/64/80 ports (on demand)
161. Electrode cap (optional):

- MRINEO 19 X
- MRINEO 31 X
- MRINEO 61 X

162. Cable for EL CFM electrodes (as required):

- CB CFM050
- CB CFM150

163. Electrodes with cable for LTM (as required):

- CB EL TPCO 19
- CB EL TPCO 31

164. Electrode (plate) EL CFM for 2-channel EEG (on demand)

165. Disk electrode with cable (as needed):

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167. Ear electrode (clip) EL TPEAR with cable (on request)

168. EL TPLINK cable for combining two ear electrodes EL TPEAR (optional)

169. Connecting cable for LTM / LTM 256 (on request)

170. CB PUD150 cable for EL PUD electrodes (optional)

171. Disposable electrodes EL PUD (on demand)

172. Monopolar needle electrode (as required):

- EMG MBOT2530
- EMG MBOT2540
- EMG MBOT3745
- EMG MBOT5050
- EMG MBOT7570

173. Snoring sensor (optional):

- EPM 1200
- EPM 1250

174. EPM 1312 breath sensor kit. includes sensor (EPM 1387) and 4 tapes in different sizes (EPM 1340 + EPM 1341 + EPM 1342 + EPM 1343) (optional)

175. Pediatric breathing sensor EPM 1390 (on demand)

176. Pediatric air flow sensor EPM 1400 (on demand)

GM / HM disposable breath sensor set , includes: EPM 1415 cable and 10 sensors adult 1426 EPMs / medium size 1427 EPMs / infant size 1428 EPMs) (as required)

178. Air flow sensor EPM 1450 (optional)

179. Set of air flow and pressure sensors EPM 1475 (optional)

180. Air pressure sensor and snoring EPM 14830 (optional)

181. Body position sensor EPM 1563 (optional)

182. EPM 1715 motion sensor kit (includes piezo transducer (EPM 1770) and pairs of ankle straps (EPMS 1730) (optional)

183. Chest force sensor kit, including interface cable and 4 chest straps (as required):

- EPM 9102 A for adults
- EPM 9102 C pediatric

184. Set of abdominal force sensor, including interface cable and 4 abdominal bands (as required):

- EPM 9101 A for adults
- EPM 9101 C pediatric

185. A set of abdominal and thoracic force sensors, including abdominal and thoracic interface cables, and 4 tapes (as required):

- EPM 914x A for adults
- EPM 914x C pediatric
- EPM 915x A for adults
- EPM 915x C pediatric

186. Tape for the breath sensor (as needed):

- EPMs 1340
- EPMs 1341
- EPMs 1342
- EPMs 1343

187. Piezocrystalline breath sensor EPMs 1387 (optional)

188. Cable for single-channel airflow sensor EPMs 1415 (optional)

189. Thermal breath sensor (on demand):

- EPMs 1426 adult size
- EPMs 1427 medium size
- EPMs 1428 pediatric

190. EPMs 1450-B interface for EPM 1450 airflow sensor (optional)
S transducer for EPM 1450 airflow sensor (optional)

192. Combined thermistor pressure and airflow sensor EPMs 1475

CO 2 cannula (as required):

- EPMs 14802, nasal and oral, male adult
- EPMs 14803, male for newborns
- EPMs 14804, male pediatric
- EPMs 14805, male adult

194. Piezoelectric pressure transducer for airflow and snoring EPMs 14830 (optional)

195. Elastic tape EPMs 1550 for body position sensor EPM 1520 (optional)

196. EPMs 1563 transducer for EPM 1563 body position sensor (optional)
197. ThermoCan cannula (as required):
 - EPMs 15802
 - EPMs 15807
198. Pair of EPMs 1730 ankle straps (optional)
199. Piezoelectric transducer EPMs 1770 for motion sensor EPM 1715 (optional)
200. A pair of induction tapes (as needed):
 - 9002-L120 EPMs
 - 9002-L133 EPMs
 - 9002-L40 EPMs
 - 9002-L60 EPMs
 - 9002-L90 EPMs
201. Adjustable induction tape EPMs 9010 (as required)
202. Induction interface cable (as required):
 - EPMs 9101 for the abdomen
 - EPMs 9102 for the thoracic region
203. Disposable ECG electrode ECG EL (optional)
204. Cable for ECG electrode (optional):
 - ECG EL TPCLB black
 - ECG EL TPCLR red

IX. ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL STUDIES MICROMED , PERFORMED BY: MORPHEUS for PSG, INCLUDED

System software PLUS EVOLUTION - BRAIN QUICK including the following software options:

Software option for registration and viewing:

- 2.1 Software option for recording and viewing EEG, 2-channel CFM version (on request)
- 2.2 Software option for recording and viewing EEG, 8-channel CFM version (on request)
- 2.3 Software option for recording and viewing EEG, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
- 2.4 Software option for recording and viewing digital video-EEG, version CLINIC (on request)
- 2.5. Software option for recording and viewing digital video-EEG in Full HD format (on request)
- 2.6 Software option for recording and viewing digital video-EEG over the network, version CLINIC (on request)
- 2.7 Software option for real-time logging and viewing over the network (as required)
- 2.8 Software option for EEG analysis in "on-line" and "off-line" modes, version: full / CFM / DSA. (on demand)
- 2.9 Software option for transmitting EEG data via TCP/IP (optional)
- 2.10 Software option " ERP -manager" (optional)
- 2.11 Software option "Management of outpatient registrar" (as required)
3. Software options for viewing:
 - 3.1. EEG software viewing option, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
 - 3.2 Software option for viewing EEG, version: EASY (on request)
 - 3.3 Software option for viewing EEG in "off-line" mode, version: full, CFM, DSA (on request)
 - 3.4 Software option for connecting to databases of external applications (optional)
 - 3.5 EEG software option "Event identification" (as required)
 - 3.6 EEG software option "Analysis of high-frequency oscillations" (on request)

- 3.7. Software option EEG "Brain Mapping", version: CLINIC, STUDIO (on request)
 - 3.8. EEG software option "Data exchange with MRI" (on request)
 - 3.9. Software EEG option "Interface HOLBERG SCORE" (optional)
 - 3.10. Software EEG option "Audio playback" (optional)
 - 3.11. Software EEG option "Averaging in "off-line" mode" (on demand)
 - 3.12. Software option for viewing digital video-EEG, version: CLINIC, STUDIO (on request)
 - 3.13. Software option for viewing digital video-EEG in FullHD format (on request)
 - 3.14. Software option for viewing digital video-EEG in multichannel mode, version: CLINIC (on request)
 - 3.15. Software option for advanced management of the archive server (optional)
 - 3.16. Program statistics option (optional)
 - 3.17. Software option for sleep program interface (optional)
 - 3.18. Software option for background data transfer (optional)
 - 3.19. Software option for exporting EEG data to the Viewer (optional)
 - 3.20. Software option to enable HIS interface (optional)
4. Software options for sleep analysis:
- 4.1 Software option for analyzing imported curves in EDF format (optional)
 - 4.2 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring SD LTM 32 EXPRESS (32 channels) (on request)
 - 4.3 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring 32 channels (on request)
 - 4.4 Software option for the analysis of curves recorded on the EEG Head Units of the SD PLUS series (on request)
5. Program control options.
- 5.1 Software option "Centralized control of video cameras" on a USB -key (on request)
 - d. Software option "Event management" on a USB -key (optional) .
6. USB key for registration and viewing programs (optional)
7. USB dongle for viewers for 1/5/10/25 users (optional)
8. USB dongle for sleep analysis programs (optional)
11. LCD display: 19"/22"/24"/27" (optional)
 12. MORPHEUS EEG head unit for PSG
 13. Electrode box JB MORPHEUS (24 channels) to EEG head unit MORPHEUS for PSG
 14. Electrode box JB MORPHEUS CLINIC (18 channels) to MORPHEUS EEG head unit for PSG (on demand)
 15. Photostimulator LED with bracket
 16. BQ NET POE M interface module for connecting MORPHEUS EEG headunit for PSG, includes SUPP BQ NET POE 2 wall mount and BQ USB interface module
 17. BQ NET POWER POE power supply for BQ NET POE series interface modules (optional)
 18. BQ USB Interface Module to MORPHEUS EEG Headunit for PSG (optional)
 19. BOX TERMINAL interface module for MORPHEUS EEG headunit for PSG, including FO PLUS circuit and SD power supply (optional)
 20. USB adapter to serial port ADAPT USB SER (optional)
 21. Isolated USB to Serial Adapter ADAPT USB SER(optional) ISO(optional)
 22. CB DC cable for connecting DC channels to MORPHEUS EEG headunit for PSG (optional)
 23. JB MORPHEUS CB cable for JB MORPHEUS electrode box
 24. CB MAG RAPID cable for connecting the Micromed system to Magstim magnetic stimulation

25. Trigger cable CB RCRC4N for stimulator (optional)
26. Trigger cable CB RCRC4R for stimulator (optional)
27. Trigger cable CB RCRC1R for stimulator (optional)
28. Pair of CB RCRC5 PAIR trigger cables for stimulator (optional)
29. CB SER RCA EXT cable for synchronization with external stimulators
30. EXTERNAL BATTERY (optional)
31. EXTERNAL CHARGER charger for external battery (on request)
32. HBC EXTERNAL BATTERY cable for connecting an external battery (optional)
33. Fiber optic cable 10m/15m/20m/30m/40m/50m/100m (optional)
34. FO ADAPTER LTM to connect BQ USB interface module to fiber optic cable (optional)
35. LED POWER power supply for LED photostimulator
36. 180 Extension cable LED EXT for LED photostimulator (on request)
37. TERM LTM Quick Disconnect Cable to EEG Head Units (optional)
38. CAMERA IP HD PRO video camera (optional)
39. CAMERA HD Full HD (1920 x 1080) camcorder with PTZ control and automatic IR switching (on demand)
40. CAMERA IP HD OUTDOOR video camera (optional)
41. CAMERA IP HD PLUS video camera (optional)
42. Camcorder CAMERA LAPTOP with USB 2.0, including microphone (optional)
43. CAMERA IP COMPACT video camera (optional)
44. GC MM video card with dual monitor output (optional)
45. HDMI EXTENDER device for transmitting HDMI signal over Cat5e/6 cable, including receiver and transmitter (optional)
46. IP MIC microphone for IP video cameras (optional)
47. IP POWER power adapter for IP video cameras (optional)
48. VLAMP IR illuminator with auto-on (on demand)
49. Composite switcher VMIXER PRO with 2 inputs (optional)
50. Professional switcher VMIXER PRO with 5 HDMI inputs (optional)
51. Bracket ARM GM MORPH to EEG head unit MORPHEUS for PSG assembly (on request)
52. Vertical bracket ARM LCD DOUBLE for two LCD displays for fixing on the cart (optional)
53. Horizontal bracket ARM LCD HORIZ for two LCD displays for fixing on the cart (optional)
54. Vertical bracket ARM LCD SW for LCD displays for fixing on the trolley (optional)
55. Adjustable BAG BELT HP to secure the headunit bag to the patient
56. BAG MORPH bag for MORPHEUS EEG head unit for PSG
57. Belt to the bag of the EEG head unit MORPHEUS for PSG
58. Rack with base BASE STAND (optional)
59. BQ CLAMP (optional)
60. Swivel mount CAM JOIN for a video camera on a stand (optional)
61. CAM STAND for mounting a video camera or head unit on a trolley (optional)
62. FIX GM mount for trolley bracket (optional)
63. Fixing FIX GM TW for table bracket (optional)
64. JOIN GM connection for attaching the head unit holder to the bracket (optional)
65. LTM SCREW screw for attaching the electrode box to the head unit (optional)
66. LED STAND PHOTIC wheel base for LED photostimulator (on demand)
67. Bracket LED SUPP PHOTIC for fixing the LED photostimulator (optional)
68. Holder STAND CAM SW video camera assembly (optional)
69. Stand STAND POLE NEW on wheels with weighting (optional)
70. SUPP BQ NET POE2 wall mount for BQ NET POE 2 / BQ NET POE M interface modules (optional)
71. SUPP MORPH Holder for MORPHEUS EEG Headunit for PSG (optional)
72. VCLAMP mount for video camera or IR light (optional)
73. VSTAND tripod stand for camcorder or IR illuminator (optional)

74. VSUPP Wall Mount for Video Camera (Optional)
75. VSUPP IP metal support for attaching an IP video camera to a swivel mount (optional)
76. VSUPP IP IR metal support for attaching the IP video camera to the rack (optional)
77. Wall mount VSUPP IP HD PLUS BLACK for video camera, black (optional)
78. VSUPP IP HD PLUS WHITE wall mount for video camera, white (optional)
79. ADAPT BNC RCA adapter (optional)
80. ADAPT DVI HDMI adapter (optional)
81. ANY WHERE 2 USB over LAN device with 2 USB ports (optional)
82. ANY WHERE 5 USB over LAN device with 5 USB ports (optional)
83. USB joystick CAM JOYSTICK for camera control (optional)
84. 24-key CAM PAD numeric keypad with USB connection for camera control (optional)
85. Numeric keypad CAM PAD 48 x 48 keys with USB connection for camera control (optional)
86. Connecting cable CB BBNB1 (BNC to BNC) (optional)
87. Connecting cable CB BNRC1 (BNC - RCA (Red) (as required)
88. Connecting cable CB BNRC5 (BNC - RCA (Black) (optional)
89. CB DVI cable for monitor (optional)
90. HDMI - HDMI MINI cable for video transmission, length 2m/5m/10m (optional)
91. HDMI to HDMI cable for high-definition video transmission, 2m/5m/15m long (optional)
92. Audio cable CB JACK JMINI with 6.3 mm plugs - MiniJack 3.5 mm (optional)
93. Audio cable CB JACK RCRC with 6.3 mm plugs - double RCA (optional)
94. Microphone cable CB MICROPHONE (optional)
95. CB MON cable for monitor (D-Sub 15 pin - D-Sub 15 pin) (optional)
96. Two-wire audio cable CB RCRC JMINI, RCA male-MiniJack 3.5 mm (optional)
97. Two-wire audio cable CB RCRC PAIR, RCA male-double RCA (as required)
98. CB VIDEO video cable as required)
99. CB IP HD PLUS 15 cable for controlling IP cameras (optional)
100. CB IP HD PLUS DC 5 cable for connecting IP HD PLUS cameras to each other (optional)
101. CB IP HD PLUS EXTENDER adapter for IP HD PLUS signal transmission over Cat5e/6 cable, pair (optional)
102. HDMI BOOSTER self-powered HDMI amplifier (optional)
103. Trolley CARSTAR mobile (on demand)
104. Bracket FIX TOR for attaching the isolation transformer to the trolley (optional)
105. Mouse pad small/regular/large (as needed)
106. Trolley CAR PLUS mobile, adjustable in height (as required)
107. FIX CART plates for fixing trolley legs (optional)
108. Power cable set KIT CB POWER for trolley (as required)
109. Set of labels KIT LABEL PLUS for the trolley (as required)
110. SUPP KEYB PLUS keyboard drawer for trolley (optional)
111. SUPP PC PLUS trolley mount (optional)
112. Small computer mount SUPP PC PLUS SFF on cart (optional)
113. Horizontal bar BAR POLE NEW for attaching the bracket to the trolley
114. Trolley CAR TORO with shelves, standard size (as required)
115. Trolley CAR TORO SMALL with shelves, small size (on demand)
116. COVER TORO rear cover for standard CARTORO trolley (optional)
117. Back cover COVER TORO SMALL for small cart CAR TORO SMALL (as required)
118. Isolation transformer TOR 500. Pmax 500 VA, black (on request)
119. Isolation transformer TOR 500 W. Pmax 500 VA, white (on request)
120. Isolation transformer TR 50, 100 VA (on request)
121. NET ISOLAT network isolator for 10/100 Mbps networks, with RJ45 connectors (optional)
122. Network cable NET RJ45 with RJ45 connector, length 3m/5m/10m (optional)
123. 4-port network switch NET SWITCH 4 (optional)
124. 8-port network switch NET SWITCH 8 (optional)

125. Professional network switch NET SWITCH PRO (optional)
126. Standard Bluetooth USB interface (optional)
127. BS USB READER for Compact Flash card reader (optional)
128. MORPH BUTTON Patient Button for use with MORPHEUS EEG Head Unit for PSG
129. PLIERS TP forceps for removing electrodes from the electrode box (optional)
130. USB extension cable USB EXT CB (optional)
131. USB EXTENDER device for transmitting a USB signal over a Cat5e / 6 cable, including a receiver and transmitter (optional)
132. BQ CONTACT contact closure control module (optional)
133. Patient button with 5m cable and RCA connector (optional)
134. BUTTON BOX to control patient buttons (optional)
135. Cables for BUTTON BOX (optional)
136. Illuminated patient button with 2 m cable and LEMO connector (optional)
137. Power supply for BUTTON BOX (optional)
138. Cable for connection unit for oximetric sensors. Length 3 m (as required)
139. Oximetry Sensor Connection Unit for use with MORPHEUS EEG Headunit for PSG (optional)
140. Oximetric sensor "EarClip", adult (on request)
141. FingerClip oximetry sensor, adult (optional)
142. Oximetric sensor "FingerClip", child (on request)
143. Fixing tape for SD OXY RA, adult size (as required)
144. Fixing tape SD OXY RC, child size (as required)
145. Flexible oximetry probe, adult size (optional)
146. Flexible oximetry sensor, child size (optional)
147. Large size reusable soft oximeter probe (digital readout 12.5 - 25.5 mm) (optional)
148. Reusable medium size soft oximeter probe (digital readout 10 - 19 mm) (optional)
149. Reusable soft oximetry sensor small size (digital display 7.5 - 12.5 mm) (optional)
150. KVM-extender COMP DVI RX compact, for receiving DVI + USB HDI (optional)
151. KVM Extender COMP DVI RX AUD compact, for receiving DVI + USB HDI + analog audio signal (optional)
152. KVM Extender COMP DVI TX compact, for DVI + USB HDI transmission (optional)
153. COMP DVI TX AUD compact KVM extender for DVI + USB HDI + analog audio transmission (optional)
154. DVI RX KVM module for DVI + USB HID reception (optional)
155. KVM Module DVI RX AUD Dual to receive DVI + USB HID + analog audio signal (as needed)
156. DVI TX KVM module for DVI + USB HID transmission (optional)
157. KVM Module DVI TX AUD Dual to Transmit DVI + USB HID + Analog Audio Signal (As Required)
158. KVM Module DVI TX AUD DCAT Dual for DVI + USB HID + Analog Audio Transmission, Dual Channel Transmission (as needed)
159. DVI TX DCAT KVM Module for DVI + USB HID Transmission, with Dual Channel Transmission (optional)
160. KIT KVM KVM module including transmitter and receiver for keyboard, mouse and monitor. (optional)
161. KVM module KIT KVM DUAL, including transmitter and receiver, for keyboard, mouse and two monitors (optional)
162. KVM Chassis KVM BODY 2 for 2 DVI modules (optional)
163. KVM Chassis KVM BODY 4 for 4 DVI modules (optional)
164. KVM BODY 4R Chassis for 4 DVI modules, including rack mount bracket (optional)
165. KVM Chassis KVM BODY 6 for 6 DVI modules (optional)
166. KVM BODY 6R Chassis for 6 DVI modules, including rack mount bracket (optional)

- 167.KVM IP Module, KVM IP RX, Receiver (optional)
- 168.KVM IP module, KVM IP TX, Transmitter (optional)
- 169.KVM matrix KVM SWITCH for 8/16/32/48/64/80 ports (on demand)
- 170. Electrode cap (as needed):

- MRINEO 19 X
- MRINEO 31 X
- MRINEO 61 X

171. Cable for EL CFM electrodes (as required):

- CB CFM050
- CB CFM150

172. Electrodes with cable for LTM (as required):

- CB EL TPCO 19
- CB EL TPCO 31

173. Electrode (plate) EL CFM for 2-channel EEG (on demand)

174. Disk electrode with cable (as required):

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176. Ear electrode (clip) EL TPEAR with cable (on request)

177.EL TPLINK cable for combining two ear electrodes EL TPEAR (on request)

178.Connecting cable for LTM / LTM 256 (on request)

179.CB PUD150 cable for EL PUD electrodes (optional)

180. Disposable electrodes EL PUD (on demand)

181. Electrode monopolar needle (on demand):

-EMG MBOT2530

- EMG MBOT3745
- EMG MBOT5050
- EMG MBOT7570
- 182. Snoring sensor (on demand):
- EPM 1200
- EPM 1250
- 183. EPM 1312 breathing sensor kit, includes sensor (EPM 1387) and 4 tapes in different sizes (EPM 1340 + EPM 1341 + EPM 1342 + EPM 1343) (optional)
- 184. Pediatric breathing sensor EPM 1390 (on request)
- 185. Pediatric airflow sensor EPM 1400 (on demand)
- GM / HM disposable breath sensor set , includes: EPM 1415 cable and 10 sensors adult 1426 EPMs / medium size 1427 EPMs / infant size 1428 EPMs) (as required)
- 187. Air flow sensor EPM 1450 (optional)
- 188. Set of air flow and pressure sensors EPM 1475 (optional)
- 189. Air pressure sensor and snoring EPM 14830 (optional)
- 190. Body position sensor EPM 1563 (optional)
- 191. EPM 1715 motion sensor kit (includes piezo transducer (EPM 1770) and pairs of ankle straps (EPMS 1730) (optional)
- 192. Chest force sensor kit, including interface cable and 4 chest straps (as required):
- EPM 9102 A for adults
- EPM 9102 C pediatric
- 193. Set of abdominal force sensor, including interface cable and 4 abdominal bands (as required):
- EPM 9101 A for adults
- EPM 9101 C pediatric
- 194. A set of abdominal and thoracic force sensors, including abdominal and thoracic interface cables, and 4 tapes (as required):
- EPM 914x A for adults
- EPM 914x C pediatric
- EPM 915x A for adults
- EPM 915x C pediatric
- 195. Tape for the breath sensor (as needed):
- EPMs 1340
- EPMs 1341
- EPMs 1342
- EPMs 1343
- 196. Piezocrystalline breath sensor EPMs 1387 (optional)
- 197. Cable for single-channel airflow sensor EPMs 1415 (optional)
- 198. Thermal breath sensor (as required):
- EPMs 1426 adult size
- EPMs 1427 medium size
- EPMs 1428 pediatric
- 199. EPMs 1450-B interface for EPM 1450 airflow sensor (optional)
- S transducer for EPM 1450 airflow sensor (optional)
- 201. Combined thermistor pressure and airflow sensor EPMs 1475
- CO 2 cannula (as required):
- EPMs 14802, nasal and oral, male adult
- EPMs 14803, male for newborns
- EPMs 14804, male pediatric
- EPMs 14805, male adult
- 203. Piezoelectric pressure transducer for airflow and snoring EPMs 14830 (optional)
- 204. Elastic tape EPMs 1550 for body position sensor EPM 1520 (optional)
- 205. EPMs 1563 transducer for EPM 1563 body position sensor (optional)
- 206. ThermoCan cannula (optional):
- EPMs 15802
- EPMs 15807

- 207. Pair of EPMs 1730 ankle straps (optional)
- 208. Piezoelectric transducer EPMs 1770 for motion sensor EPM 1715 (optional)
- 209. A pair of induction tapes (as needed):
 - 9002-L120 EPMs
 - 9002-L133 EPMs
 - 9002-L40 EPMs
 - 9002-L60 EPMs
 - 9002-L90 EPMs
- 210. Adjustable induction band EPMs 9010 (optional)
- 211. Induction interface cable (optional):
 - EPMs 9101 for the abdomen
 - EPMs 9102 for the thoracic region
- 212. Inductive interface cable with connector for MORPHEUS (optional):
 - EPMs 914x
 - EPMs 915x
- 213. Disposable ECG electrode ECG EL (on request)
- 214. Cable for ECG electrode (optional):
 - ECG EL TPCLB black
 - ECG EL TPCLR red

1.2. Main parameters and characteristics

1.2.1 Classification and description of the EEG system

1.2.1. 1 The EEG system consists of electronic computing and electrical equipment, blocks, amplifiers, electrodes, brackets, racks, cables, accessories and software pre-installed in the system (the list of all products that form the system, depending on the version of the system, is presented in the **Appendix A (Included in the Package)**), suitable for working with the *Microsoft operating system Windows XP* and older. The system versions differ ONLY in the number of channels in the EEG blocks, all components of the systems are identical and have the same technical and functional characteristics. In general, the system consists of:

System software PLUS EVOLUTION - BRAIN QUICK including the following software options:

Software option for registration and viewing:

- 2.1 Software option for recording and viewing EEG, 2-channel CFM version (on request)
- 2.2 Software option for recording and viewing EEG, 8-channel CFM version (on request)
- 2.3 Software option for recording and viewing EEG, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
- 2.4 Software option for recording and viewing digital video-EEG up to 64 channels, version CLINI C (on request)
- 2.5 Software option for recording and viewing digital video-EEG, version CLINI C (on request)
- 2.6. Software option for recording and viewing digital video-EEG in Full HD format (on request)
- 2.7 Software option for recording and viewing digital video-EEG over the network, version CLINI C (on request)
- 2.8 Software option for recording and viewing digital video-EEG over a network up to 64 channels, version CLINI C (on request)
- 2.9 Software option for logging and viewing over the network in real time (as required)
- 2.10 Software option for EEG analysis in "on-line" and "off-line" modes, version: full / CFM /

- 3.1. EEG software viewing option, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
- 3.2 Software option for viewing EEG, version: EASY (on request)
- 3.3 Software option for viewing EEG in "off-line" mode, version: full, CFM, DSA (on request)
- 3.4 Software option for connecting to databases of external applications (optional)
- 3.5 EEG software option "Event identification" (as required)
- 3.6 EEG software option "Analysis of high-frequency oscillations" (on request)
- 3.7. Software option EEG "Brain Mapping", version: CLINIC, STUDIO (on request)
- 3.8. EEG software option "Data exchange with MRI" (on request)
- 3.9. Software EEG option "Interface HILBERG SCORE" (optional)
- 3.10. Software EEG option "Audio playback" (optional)
- 3.11. Software EEG option "Averaging in "off-line" mode" (on demand)
- 3.12. Software option for viewing digital video-EEG, version: CLINIC, STUDIO (on request)
- 3.13. Software option for viewing digital video-EEG in Full HD format (on request)
- 3.14. Software option for viewing digital video-EEG in multichannel mode, version: CLINIC (on request)
- 3.15. Software option for advanced management of the archive server (optional)
- 3.16. Program statistics option (optional)
- 3.17. Software option for sleep program interface (optional)
- 3.18. Software option for background data transfer (optional)
- 3.19. Software option for exporting EEG data to the Viewer (optional)
- 3.20. Software option to enable HIS interface (optional)

- 4.1 Software option for analyzing imported curves in EDF format (optional)
- 4.2 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring SD LTM 32 EXPRESS (32 channels) (on request)
- 4.3 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring SD LTM 64 EXPRESS (64 channels) (optional)
- 4.4 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring 32 channels (on demand)
- 4.5 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring 64 channels (on request)
- 4.6 Software option for the analysis of curves recorded on the EEG Head Units of the SD PLUS series (on request)

5.1 Software option "Centralized control of video cameras" on a USB -key (on request)
5.2. Software option "Event management" (on demand) on a USB key

8. USB dongle for sleep analysis programs (optional)

80

9 CB TERMINAL 15 cable for connecting the interface board to the BOX TERMINAL EXPRESS interface module

EXPRESS wall panel WALL CONNECTOR /12 with additional socket and IN-trigger and OUT-trigger for use in combination with BOX TERMINAL EXPRESS interface module (optional) and 12V voltage amplifier (optional)

11LCD: 19"/22"/24"/27" (optional)

12EEG head unit for long-term monitoring SD LTM 32 EXPRESS (32 channels)

13 Electrode box JB LTM 32 EXPRESS to the EEG head unit for long-term monitoring SD LTM 32 EXPRESS (32 channels)

14BQ interface module USB PLUS to EEG headunit for long-term SD monitoring LTM PLUS (64 channels)

15EEG head unit for long-term SD monitoring LTM 32 PLUS (32 channels)

16JB LTM 32 PLUS Electrode box to the EEG head unit for long-term SD monitoring LTM 32 PLUS (32 channels)

17JBD electrode box SUB PLUS with Dual 37-Pin D - Sub Connector to EEG Head Units for SD Series Long-Term Monitoring LTM PLUS

18 JB TOPHAT PLUS Electrode Box for all channels with dual 37-pin D-Sub connector to EEG Head Units for long-term monitoring SD series LTM PLUS

19 BOX TERMINAL MULTI interface module for connecting the EEG head unit for long-term SD monitoring LTM PLUS (up to 256 channels)

20BQ NET POE PLUS interface module for connecting EEG head unit for SD long-term monitoring LTM PLUS , includes SUPP wall mount BQ NET POE 2

21 CB LTM DC cable for connecting DC channels to EEG Head Units for long-term SD monitoring LTM

22EEG head unit for long-term monitoring SD LTM 64 EXPRESS (64 channels)

23JB LTM 64 EXPRESS electrode box to the EEG head unit for long-term monitoring SD LTM 64 EXPRESS (64 channels)

SD EEG Head Units for Long-Term Monitoring LTM (64 channels)

25 Electrode box JB LTM 64 PLUS to the EEG head unit for long-term monitoring SD LTM 64 PLUS (64 channels)

26EEG head unit for long-term monitoring SD LTM 64 PLUS (64 channels)

27EEG head unit SD PLUS CLINIC (29 channels)

28EEG head unit SD PLUS FLEXI (38 channels)

29EEG head unit SD PLUS HIGH RATE (29 channels)

30EEG head unit SD PLUS RESEARCH (48 channels)

31 MORPHEUS EEG head unit for PSG

32 JB MORPHEUS Electrode Box (24 channels) to MORPHEUS EEG Head Unit for PSG

33BQ NET POE M interface module for connecting MORPHEUS EEG headunit for PSG. includes SUPP BQ NET POE 2 wall mount and BQ USB interface module

34BQ USB Interface Module to MORPHEUS EEG Headunit for PSG (as required)

- 35 BOX TERMINAL interface module for MORPHEUS EEG head unit for PSG, including FO PLUS circuit and SD power supply (on request)
- 36 CB DC cable for connecting DC channels to MORPHEUS EEG headunit for PSG (optional)
- 37 CB JB MORPHEUS cable for JB MORPHEUS electrode box
- 38 ARM GM MORPH bracket for MORPHEUS EEG head unit for PSG assembly (on request)
- 39 BAG MORPH bag for MORPHEUS EEG head unit for PSG
- 40 Belt to the bag of MORPHEUS EEG head unit for PSG
- 41 SUPP MORPH Holder for MORPHEUS EEG Head Unit for PSG (optional)
- 42 MORPH BUTTON Patient Button for use with the 43 MORPHEUS EEG Headunit for PSG
- 44 Oximetry sensor connection box for use with MORPHEUS EEG headunit for PSG (optional)
- 45 JB LTM DSUB Electrode Box with Dual 37-Pin D-Sub Connector for EEG Head Units for Long-Term Monitoring SD LTM EXPRESS Series (optional)
- 46 Electrode box JB LTM TOPHAT (universal) with sensor inputs for all channels (as required)
- 47 LED photostimulator with bracket
- 48 BQ NET POE 2 interface module for connecting SD LTM EXPRESS long-term monitoring EEG headunit, includes SUPP BQ NET POE 2 wall mount and BQ USB EXPRESS interface module (optional)
- EXPRESS wall panel WALL CONNECTOR / DE with additional socket and IN-trigger and OUT-trigger for use in combination with BOX TERMINAL EXPRESS interface module (optional)
- 50 BQ NET POWER POE power supply for BQ NET POE series interface modules (optional)
- 51 BQ PCI EXPR LTM FULL interface board based on PCI Express FULL Size to control EEG headunit for long-term monitoring SD LTM EXPRESS (optional)
- 52 BQ PCI EXPR LTM HALF interface board based on PCI Express HALF Size to control the EEG headunit for long-term monitoring SD LTM EXPRESS (on request)
- 53 BQ USB EXPRESS interface module to EEG head unit for long-term monitoring SD LTM EXPRESS (on request)
- 54 BQ USB EXPRESS interface module to EEG head unit for long-term monitoring SD LTM EXPRESS (optional)
- 55 BQ USB PLUS interface module to EEG head unit for long-term monitoring SD LTM EXPRESS (on request)
- 56 USB adapter to serial port ADAPT USB SER (optional)
- 57 Isolated USB to Serial Adapter ADAPT USB SER ISO (as needed)
- 58 CB LTM DC cable for connecting DC channels to EEG Head Units for long-term SD monitoring LTM
- 59 CB MAG RAPID cable for connecting the Micromed system to **Magstim magnetic stimulation devices** (D-Sub 26 pin (optional)
- 60 Trigger cable CB RCRC4N for stimulator (optional)
- 61 Trigger cable CB RCRC4R for stimulator (optional)
- 62 CB RCRC1R trigger cable for stimulator (optional)

- 63 Pair of CB RCRC5 PAIR trigger cables for stimulator (optional)
- 64 CB SER RCA EXT cable for synchronization with external stimulators (optional)
- 65 Fiber optic cable 10m/15m/20m/30m/40m/50m/100m (optional)
- 66 HBC EXPRESS cable for connecting an interface board or an interface module to EEG Head Units (as required)
- 67 Cable HBC INT LTM R for connecting an interface board or interface module to EEG Head Units for long-term monitoring of the SD LTM series (on request)
- 68 LED POWER power supply for LED photostimulator
- 69 TERM LTM Quick Disconnect Cable to EEG Head Units
- 70CAMERA IP HD PRO (optional)
- 71 CAMERA HD Full HD (1920 x 1080) with PTZ control and automatic IR switching (on demand) (on demand)
- 72 CAMERA IP HD OUTDOOR (optional)
- 73 CAMERA IP HD PLUS (optional)
- 74 CAMERA LAPTOP with USB 2.0, including microphone (optional)
- 75 CAMERA IP COMPACT video camera (optional)
- 76 GC MM video card with dual monitor outputs (optional)
- 77 HDMI EXTENDER device for transmitting HDMI signal over Cat 5 e /6 cable, including receiver and transmitter (optional)
- 78 IP MIC for IP cameras (optional)
- 79 IP POWER power adapter for IP cameras (optional)
- 80 VLAMP IR auto-on (on demand)
- 81 Composite switcher VMIXER PRO with 2 inputs (optional)
- 82 VMIXER PRO professional switcher with 5 HDMI inputs (optional)
- 83 Bracket ARM GM LTM for EEG head units for long-term monitoring of the SD LTM series, complete (on request)
- 84 Vertical bracket ARM LCD DOUBLE for two LCD displays for fixing on the cart (optional)
- 85Horizontal bracket ARM LCD HORIZ for two LCD displays for fixing on the cart (optional)
- 86 ARM LCD SW Vertical LCD Bracket for Trolley Fixing (Optional)
- 87 Adjustable BAG BELT HP for fixing the headunit bag to the patient
- 88 EXPRESS WALL CONNECTOR NT /12 wall panel with additional socket for use in combination with BOX TERMINAL interface module (optional)EXPRESS and 12V voltage booster (optional)
- 89 BAG LTM HP Bag for EEG Head Unit for Long-Term Monitoring SD LTM Series (optional)
- 90 BASE STAND (optional)
- 91 BQ CLAMP (optional)
- 92 CAM JOIN swivel mount for video camera on a stand (optional)

- 93 CAM STAND for mounting a video camera or head unit on a trolley (optional)
- 94 FIX GM attachment for trolley bracket (optional)
- 95 Fixture FIX GM TW for table bracket (optional)
- 96 JOIN GM connection for attaching the head unit holder to the bracket (optional)
- 97 LTM SCREW screw for attaching the electrode box to the head unit (optional)
- 98 LED STAND PHOTIC wheel base for LED photostimulator (on request)
- 99 LED SUPP PHOTIC bracket for fixing the LED photostimulator (optional)
- 100 STAND CAM SW camera holder assembly (optional)
- 101 STAND POLE NEW on wheels with weight (optional)
- 71 Wall mount SUPP BQ NET POE 2 wall mount for BQ NET POE 2 interface modules BQ NET POE 2 / BQ NET POE M (optional) (optional) (optional)
- 102 SUPP LTM Holder for SD LTM Series EEG Head Units (optional)
- 103 VCLAMP mount for video camera or IR light (optional)
- 104 VSTAND tripod stand for video camera or IR light (optional)
- 105 VSUPP wall bracket for video camera (optional)
- 106 VSUPP IP metal support for attaching the IP video camera to the swivel mount (optional)
- 107VSUPP IP IR metal support for attaching the IP video camera to the rack (optional)
- 108 VSUPP IP HD PLUS BLACK wall mount for video camera, black (optional)
- 109 VSUPP IP HD PLUS WHITE wall mount for video camera, white (optional)
- 110 ADAPT BNC RCA adapter (optional)
- 111 ADAPT DVI HDMI adapter (optional)
- 112 ANY WHERE 2 USB over LAN device with 2 USB ports (optional) (optional)
- 113 ANY WHERE 5 USB over LAN device with 5 USB ports (optional)
- 114 USB CAM JOY STICK to control cameras (optional)
- 115 24-key CAM PAD numeric keypad with USB connection for camera control (optional)
- 116CAM PAD Numeric Keypad 48 x 48 keys with USB connection to control cameras (optional)
- 117Connecting cable CB BBNB1 (BNC - BNC) (optional)
- 118 Connecting cable CB BNRC1 (BNC - RCA (Red) (optional)
- 119 Connecting cable CB BNRC5 (BNC - RCA (Black) (optional)
- 120CB DVI monitor cable (optional)
- 121HDMI - HDMI MINI cable for video transmission, length 2m/5m/10m (optional)
- 122 HDMI to HDMI cable for high-definition video transmission, 2m/5m/15m (optional)
- 123Audio cable CB JACK JMINI with 6.3 mm plugs - MiniJack 3.5 mm (optional)
- 124Audio cable CB JACK RCRC with 6.3 mm plugs - double RCA (optional)

- 125 CB MICROPHONE cable (optional)
- 126 CB MON cable for monitor (D-Sub 15 pin - D-Sub 15 pin) (optional)
- 127 Two-wire audio cable CB RCRC JMINI, RCA male to MiniJack 3.5 mm (optional)
- 128 Two-wire audio cable CB RCRC PAIR, RCA male-double RCA (as required)
- 129 CB VIDEO video cable as required)
- 130CB IP HD PLUS 15 cable for controlling IP HD PLUS cameras (optional)
- 131 CB IP HD PLUS DC 5 cable for connecting IP HD PLUS cameras to each other (optional)
- 132 CB IP HD PLUS EXTENDER adapter for IP HD PLUS signal transmission over Cat5e/6 cable, pair (optional)
- 133 HDMI BOOSTER self-powered HDMI amplifier (optional)
- 134 CARSTAR mobile trolley (on request)
- 135 FIX TOR bracket for attaching the isolation transformer to the trolley (optional)
- 136 Mouse pad small/regular/large (optional)
- 137 Trolley CAR PLUS mobile, adjustable in height (as required)
- 138 FIX CART plates for fixing trolley legs (optional)
- 139 Power cable kit KIT CB POWER for trolley (optional)
- 140 Set of labels KIT LABEL PLUS for trolley (as required)
- 141 SUPP KEYB PLUS keyboard drawer for trolley (optional)
- 142 SUPP PC PLUS trolley mount (optional)
- 143 SUPP PC PLUS SFF small computer trolley mount (optional)
- 144 Horizontal bar BAR POLE NEW for attaching the bracket to the trolley (optional)
- 145 CAR TORO trolley with shelves, standard size (as required)
- 146CAR TORO SMALL trolley with shelves, small size (on request)
- 147 COVER TORO rear cover for standard CARTORO trolley (optional)
- 148 Rear cover COVER TORO SMALL for small trolley CAR TORO SMALL (optional)
- 149 Isolation transformer TOR 500, Pmax 500 VA, black (on request)
- 150 Isolation transformer TOR 500 W, Pmax 500 VA, white (on request)
- 151Isolation transformer TR 50, 100 VA (optional)
- 152 NET ISOLAT network isolator for 10/100 Mb networks, with RJ45 connectors (optional)
- 153 Network cable NET RJ45 with RJ45 connector, length 3m/5m/10m (as required)
- 154 NET SWITCH 4 4-port network switch (optional)
- 155 NET SWITCH 8 8-port network switch (optional)
- 156 NET SWITCH PRO professional network switch (optional)
- 157 Standard Bluetooth USB interface (optional)

- 158 BS USB READER Compact Flash card reader (optional)
- 159 Extension LTM BUTTON EXT for patient button LTM BUTTON TP, length 120 cm (optional)
- LTM BUTTON TP Patient Button for use with EEG Head Units for Long-Term Monitoring SD LTM Series
- 161 PLIERS TP forceps for removing electrodes from the electrode box (optional)
- 162 USB extension cable USB EXT CB (optional)
- 163 USB EXTENDER device for transmitting USB signal over Cat5e/6 cable, including receiver and transmitter (optional)
- 164 BQ CONTACT contact closure control module (optional)
- 165 Patient button with 5 m cable and RCA connector (optional)
- 166 BUTTON BOX for controlling patient buttons (optional)
- 167 Cables for BUTTON BOX (optional)
- 168 Illuminated patient button with 2 m cable and LEMO connector (optional)
- 169 Power supply for BUTTON BOX (optional)
- 170 Cable for connection unit for oximetry sensors. Length 3 m (as required)
- 171 Oximetry Sensor Connection Box for use with EEG Head Units for long-term monitoring SD LTM series (as required)
- 172 Oximetry sensor "EarClip", adult (optional)
- 173 FingerClip oximetry sensor, adult (optional)
- 174 FingerClip oximetry sensor, child (optional)
- 175 Fixing tape SD OXY RA, adult size (on request)
- 176 Fixing tape SD OXY RC, child size (as required)
- 177 Flexible oximetry probe, adult size (optional)
- 178 Flexible oximetry sensor, child size (optional)
- 179 Large reusable soft oximeter probe (digital readout 12.5 - 25.5 mm) (optional)
- 180 Medium size reusable soft oximeter probe (digital readout 10 - 19 mm) (optional)
- 181 Reusable soft oximeter probe small size (digital readout 7.5 - 12.5 mm) (optional)
- 182 COMP DVI RX KVM extender compact, for DVI + USB HDI reception (optional)
- COMP DVI RX AUD 183KVM extender compact, for receiving DVI + USB HDI + analog audio signal (optional)
- 184 KVM Extender COMP DVI TX compact, for DVI + USB HDI transmission (optional)
- 185 KVM Extender COMP DVI TX AUD compact, for DVI + USB HDI + analog audio transmission (optional)
- 186 KVM DVI module RX for DVI + USB reception HID (on demand)
- 187 KVM DVI module RX AUD dual for DVI + USB reception HID + analog audio signal (optional)
- 188 KVM DVI module TX for DVI + USB transmission HID (on demand)

189 KVM DVI module TX AUD dual for DVI + USB transmission HID + analog audio signal (optional)

190 KVM DVI module TX AUD DCAT dual for DVI + USB transmission HID + analog audio signal with two-channel transmission (as required)

191 KVM DVI module TX DCAT to DVI + USB transmission HID ,dual transmission (optional)

192 KIT KVM module including transmitter and receiver for keyboard, mouse and monitor. (optional)

193KVM KIT KVM DUAL module, including transmitter and receiver, for keyboard, mouse and dual monitors (optional)

194 KVM -KVM Chassis BODY 2 for 2 DVI modules (optional)

195 KVM - KVM Chassis BODY 4 for 4 DVI modules (optional)

196 KVM -KVM Chassis BODY 4 R for 4 DVI modules , including rack mounting bracket (optional)

197 KVM -KVM Chassis BODY 6 for 6 DVI modules (optional)

198 KVM -KVM Chassis BODY 6 R for 6 DVI modules , including rack mounting bracket (optional)

199 IP KVM Module, KVM IP RX, Receiver (optional)

200 KVM IP, KVM IP TX, Transmitter (optional)

201 KVM matrix KVM SWITCH for 8/16/32/48/64/80 ports (on demand)

202Additional SUPP TOPHAT mount for JB LTM TOPHAT Electrode box (on request)

203 Wall panel EXPRESS WALL CONNECTOR NT /DE with additional socket for use in combination with BOX TERMINAL EXPRESS interface module (on request)

204 HBC BT/RJ45 extension adapter for BOX TERMINAL EXPRESS interface module (module side) (optional)

205 HBC HT/RJ45 extension adapter for BOX TERMINAL EXPRESS interface module (PC side) (optional)

206 EXTERNAL BATTERY (optional)

207 EXTERNAL CHARGER for external battery (optional)

208 HBC EXTERNAL BATTERY cable for external battery connection (optional)

209 Extension cable LED EXT for LED photostimulator (optional)

210 TERM EXPRESS quick-disconnect cable to the EEG Head Unit for long-term monitoring SD LTM EXPRESS in a configuration of up to 64 channels (on request)

211 Electrode cap (optional):

- MRINEO 19 X

- MRINEO 31 X

- MRINEO 61 X

212. Cable for EL CFM electrodes (as required):

- CB CFM050

- CB CFM150

213. Electrodes with cable for LTM (as required):

- CB EL TPCO 19

- CB EL TPCO 31

214. Electrode (plate) EL CFM for 2-channel EEG (on demand)

215. Disk electrode with cable (as needed):

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217. Ear electrode (clip) EL TPEAR with cable (optional)

218. EL TPLINK cable for combining two ear electrodes EL TPEAR (on request)

219. Connecting cable for LTM / LTM 256 (on request)

220. CB PUD150 cable for EL PUD electrodes (optional)

221. Disposable electrodes EL PUD (on demand)

222. Electrode monopolar needle (on demand):

- EMG MBOT2530
- EMG MBOT2540
- EMG MBOT3745
- EMG MBOT5050
- EMG MBOT7570

223. Snoring sensor (optional):

- EPM 1200
- EPM 1250

224. EPM 1312 respiration sensor set, includes sensor (EPM 1387) and 4 tapes of various sizes (EPM 1340 + EPM 1341 + EPM 1342 + EPM 1343) (as required)

225. Pediatric breath sensor EPM 1390 (on request)

226. Pediatric air flow sensor EPM 1400 (on request)

GM / HM disposable breath sensor kit , includes: EPM 1415 cable and 10 sensors adult 1426 EPMs / medium size 1427 EPMs / infant size 1428 EPMs) (as required)

228. EPM 1450 air flow sensor (optional)

229. Set of air flow and pressure sensors EPM 1475 (as required)

sensor EPM 14830 (optional)

231. Body position sensor EPM 1563 (optional)
232. EPM 1715 motion sensor kit (including piezoelectric transducer (EPM 1770) and pairs of ankle straps (EPMS 1730) (optional)
233. Chest force sensor kit, including interface cable and 4 chest straps (as required):
- EPM 9102 A for adults
 - EPM 9102 C pediatric
234. Abdominal force sensor kit, including interface cable and 4 abdominal bands (as required):
- EPM 9101 A for adults
 - EPM 9101 C pediatric
235. A set of abdominal and thoracic force sensors, including abdominal and thoracic interface cables, and 4 tapes (as required):
- EPM 914x A for adults
 - EPM 914x C pediatric
 - EPM 915x A for adults
 - EPM 915x C pediatric
236. Tape for the breath sensor (as needed):
- EPMs 1340
 - EPMs 1341
 - EPMs 1342
 - EPMs 1343
237. Piezocrystalline breath sensor EPMs 1387 (optional)
238. Cable for single-channel airflow sensor EPMs 1415 (optional)
239. Thermal breath sensor (as required):
- EPMs 1426 adult size
 - EPMs 1427 medium size
 - EPMs 1428 pediatric
240. EPMS 1450-B interface for EPM 1450 air flow sensor (optional)
- S transducer for EPM 1450 airflow sensor (optional)
242. Combined thermistor pressure and air flow sensor EPMs 1475
- CO 2 cannula (as required):
- EPMs 14802, nasal and oral, male adult
 - EPMs 14803, male for newborns
 - EPMs 14804, male pediatric
 - EPMs 14805, male adult
244. Piezoelectric pressure transducer for airflow and snoring EPMs 14830 (optional)
245. Elastic tape EPMs 1550 for body position sensor EPM 1520 (optional)
246. EPMs 1563 transducer for EPM 1563 body position sensor (optional)
247. ThermoCan cannula (optional):
- EPMs 15802
 - EPMs 15807
248. Pair of EPMs 1730 ankle straps (optional)
249. Piezoelectric transducer EPMs 1770 for motion sensor EPM 1715 (optional)
250. A pair of induction tapes (as needed):
- 9002-L120 EPMs
 - 9002-L133 EPMs
 - 9002-L40 EPMs
 - 9002-L60 EPMs
 - 9002-L90 EPMs
251. Adjustable induction band EPMs 9010 (optional)
252. Induction interface cable (as required):
- EPMs 9101 for the abdomen
 - EPMs 9102 for the thoracic region
253. Disposable ECG electrode ECG EL (on demand)
254. Cable for ECG electrode (optional):
- ECG EL TPCLB black
 - ECG EL TPCIR red

! The production of system versions is carried out according to one technical documentation, according to one technological and production process, the production quality management system is implemented and controlled by ONE for all system versions.!

1.2.1.2 The design of the EEG system and its components must comply with the working drawings.

1.2.1.3 The design of the EEG system ensures continuous operation.

1.2.1.4 The EEG system refers to active diagnostic non-invasive medical devices for long-term use, the degree of potential risk of use is 2b.

1.2.1.5 The EEG system is multi-channel computerized .

1.2.2 The EEG system provides the removal of diagnostic parameters, their transfer to a personal computer, storage and display on a liquid crystal display. The presence of oximetric sensors in the system does not classify the system or its parts as measuring instruments, since the purpose of the product is the graphic recording of bioelectric signals. The data obtained during the operation of the sensors serve for control, and do not carry a measuring function, in addition, the presence of sensors is a secondary additional functionality of the product, the product does not include an oximeter.

Leads attached to the scalp and earlobes transmit electrical signals that are further converted into an electroencephalogram .

1.2.3 The general operational and technical characteristics of the EEG system must comply with Table 1. Data on the operational and technical characteristics, depending on the version of the EEG unit, are indicated in Table 1A.

For all technical parameters specified in the text of the document (including Annexes), where not otherwise indicated, the tolerance is +/-10%.(Normal)

T a b l e 1. General technical characteristics of the complete system assembly.

Parameter name	Norm
1	2
Rated supply voltage for electronic computing equipment, V - with tolerance, %	220 (single phase) <u>+10</u>
Mains AC frequency, Hz	50
Maximum consumption current, A	16
Power consumption, kW, no more	5

Table 1 continued

1	2
Input voltage range, mV	0.03...5
Weight, kg, no more	according to working drawings
Dimensions (length, width, height)*	
Time of readiness for work after switching on (establishment of the operating mode), min., no more	5
The operating mode of the system (any version) is long. Time, h, no more	8
Input impedance, MOhm, not less	5
Input signal level, mV	0.5...5
Coefficient of attenuation of common-mode signals, not less than	30,000
Skin-galvanic EMF of electrodes, mV, not more than	<u>+ 300</u>
Common-mode signal on the electrodes relative to the total, V, not more than	1.5
Maximum value of probing current, mA	2
Maximum value of probing voltage with connected electrodes, not more than, V	1
Direct current in the patient's circuit, flowing through any electrode, excluding neutral, μ A, not more than	0.1
Identity of EEG Signal Shapes in Range and Relative Error in Measuring Amplitude and Time Parameters of EEG Signals	ensured
The level of internal noise, reduced to the input, μ V, no more	20
Relative voltage measurement error - in the range from 0.1 to 0.5 mV, %, no more; - in the range from 0.5 to 4 mV, %, no more	<u>+15</u> <u>+ 7</u>
Nonlinearity, %	<u>+ 2</u>
Sensitivity, mm/mV - relative error of sensitivity setting, %	in row 1, 2, 5, 10, 20 <u>+5</u>
Time constant, s, not less	3.2
Unevenness of the amplitude-frequency characteristic (AFC), %, in the frequency range - from 0.5 to 60 Hz; - 60 to 75 Hz	from minus 10 to plus 5

	from minus 30 to plus 5
Possible number of measurements per hour	not regulated
Relative error of measurement of time intervals, %, no more	± 7
Degree of protection against liquid penetration (unless otherwise specified in this document)	1 PX0

Table 1A.

SD Plus RESEARCH

-40 common reference channels (channels 1 - 40). 16 of them can be switched to bipolar mode (channels 25 - 40).

-02 bipolar channels, which can also be used to receive high-frequency PE (c. EP1, EP2).

-02 channels of inductive tapes (Thor, Abdo).

-03 channels for block for oximetric sensors (plethysmogram, heart rate).

-01 marker channel

- EEG channels:

- Sampling frequency varies from 256 Hz to 16 kHz per channel.

- Maximum bandwidth DC-F, where F depends on the sampling rate F_s ($0.27 \cdot F_s$).

- Full range: $\pm 3.2\text{mV}$ to $\pm 1638.4\text{mV}$.

- ADC 24 bits (resolution up to 0.38 nV).

- VP channels:

- Sampling frequency varies from 256Hz to 128KHz per channel.

- Maximum bandwidth 0.11 Hz - 12.8 kHz; the upper limit of the frequency range drops to $0.31 \cdot F_s$ if the sampling frequency drops below 64 kHz.

- Full range: $\pm 0.8\text{mV}$ to $\pm 51.2\text{mV}$.

- ADC 24 bits (resolution up to 0.095 nV).

- Gain and filtering are set using the software.

- Touch-proof connectors according to DIN 42802 for electrodes and 25-pin D-sub connector for pre-mounted headphones.

- Impedance control on each input jack, indicated by LEDs.

- 1 serial port for oximetry sensors (3 channels: plethysmogram, heart rate, does not require an external power supply).

- Internal microcontroller.

- Graphic display for displaying information and warnings (modified options are possible, limited by the software in the number of channels or performance).

- Power supply from PC via interface; isolation of up to 2500 volts AC with electromagnetic couplings on communication lines and up to 3 kV DC with a DC/DC converter on power lines, for patient safety.

SD Plus Flexi

-32 common reference channels (ch. 1 - 32) (modified options are possible, limited by the software in the number of channels or performance). 8 of them can be switched to bipolar mode (k. 25 - 32).

-02 bipolar channels, which can also be used to receive high-frequency PE (c. EP1, EP2).

-03 channels for block for oximetric sensors (plethysmogram, heart rate).

-01 marker channel

- EEG channels:

- Sampling frequency varies from 256 Hz to 16 kHz per channel.

- Maximum bandwidth DC-F, where F depends on the sampling rate F_s ($0.27 \cdot F_s$).

- Full range: $\pm 3.2\text{mV}$ to $\pm 1638.4\text{mV}$.

- ADC 24 bits (resolution up to 0.38 nV).

- VP channels:

- Sampling frequency varies from 256Hz to 128KHz per channel.

- Maximum bandwidth 0.11 Hz - 12.8 kHz; the upper limit of the frequency range drops to $0.31 \cdot F_s$ if the sampling frequency drops below 64 kHz.

- Full range: $\pm 0.8\text{mV}$ to $\pm 51.2\text{mV}$.

- ADC 24 bits (resolution up to 0.095 nV).

- Gain and filtering are set using the software.

- Touch-proof connectors according to DIN 42802 for electrodes and 25-pin D-sub connector for pre-mounted headphones.

- Impedance control on each input jack, indicated by LEDs.

- 1 serial port for the block for oximetry sensors (3 channels: plethysmogram, heart rate, does not require an external power supply).

- Internal microcontroller.

- Graphic display for displaying information and warnings (modified options are possible, limited by the software in the number of channels or performance).

- Power supply from PC via interface; isolation of up to 2500 volts AC with electromagnetic couplings on communication lines and up to 3 kV DC with a DC/DC converter on power lines, for patient safety.

SD Plus Flexi HighRate

- 29 reception channels (they can be selected on the settings screen of the reception program)

- variable sampling rate 256 Hz - 1 kHz per channel

- EEG channels:

- Sampling frequency varies from 256 Hz to 16 kHz per channel.

- Maximum bandwidth DC-F, where F depends on the sampling rate F_s ($0.27 \cdot F_s$).

- Full range: $\pm 3.2\text{mV}$ to $\pm 1638.4\text{mV}$.

- ADC 24 bits (resolution up to 0.38 nV).

- VP channels:

- Sampling frequency varies from 256Hz to 128KHz per channel.

- Maximum bandwidth 0.11 Hz - 12.8 kHz; the upper limit of the frequency range drops to $0.31 \cdot F_s$ if the sampling frequency drops below 64 kHz.

- Full range: $\pm 0.8 \text{ mV}$ to $\pm 51.2 \text{ mV}$.

- ADC 24 bits (resolution up to 0.095 nV).

- Gain and filtering are set using the software.
- Touch-proof connectors according to DIN 42802 for electrodes and 25-pin D-sub connector for pre-mounted headphones.
- Impedance control on each input jack, indicated by LEDs.
- 1 serial port for the block for oximetry sensors (3 channels: plethysmogram, heart rate, does not require an external power supply).
- Internal microcontroller.
- Graphic display for displaying information and warnings (modified options are possible, limited by the software in the number of channels or performance).
- Power supply from PC via interface; isolation of up to 2500 volts AC with electromagnetic couplings on communication lines and up to 3 kV DC with a DC/DC converter on power lines, for patient safety.

SD Plus Flexi Clinic

- 29 reception channels (they can be selected on the settings screen of the reception program)

• variable sampling rate 256 Hz - 1 kHz per channel

• EEG channels:

- Sampling frequency varies from 256 Hz to 16 kHz per channel.

- Maximum bandwidth DC-F, where F depends on the sampling rate F_s ($0.27 \cdot F_s$).

- Full range: $\pm 3.2 \text{ mV}$ to $\pm 1638.4 \text{ mV}$.

- ADC 24 bits (resolution up to 0.38 nV).

• VP channels:

- Sampling frequency varies from 256 Hz to 128 kHz per channel.

- Maximum bandwidth 0.11 Hz - 12.8 kHz; the upper limit of the frequency range drops to $0.31 \cdot F_s$ if the sampling frequency drops below 64 kHz.

- Full range: $\pm 0.8 \text{ mV}$ to $\pm 51.2 \text{ mV}$.

- ADC 24 bits (resolution up to 0.095 nV).

- Gain and filtering are set using the software.
- Touch-proof connectors according to DIN 42802 for electrodes and 25-pin D-sub connector for pre-mounted headphones.
- Impedance control on each input jack, indicated by LEDs.
- 1 serial port for the block for oximetry sensors (3 channels: plethysmogram, heart rate, does not require an external power supply).
- Internal microcontroller.
- Graphic display for displaying information and warnings (modified options are possible, limited by the software in the number of channels or performance).

- Power supply from PC via interface; isolation of up to 2500 volts AC with electromagnetic couplings on communication lines and up to 3 kV DC with a DC/DC converter on power lines, for patient safety.

SD PLUS RESEARCH, SD PLUS FLEXI, SD PLUS FLEXI HIGH RATE, SD PLUS FLEXI CLINIC:

ALL ME CLASS II APPLIED PART TYPE BF according to EN 60601-1.

To ensure compliance, use with a PC connected to the mains through an isolating or isolating transformer conforming to IEC 60742 or IEC 60989 and EN 60601-1, or with a laptop powered by a battery.

The device is intended for continuous use. The device is not suitable for use in the presence of flammable anesthetic mixtures with air, oxygen or nitrous oxide.

CHARACTERISTICS OF CHANNELS:

SD PLUS RESEARCH - 48 channels total

SD PLUS FLEXI - 38 channels Total

SD PLUS FLEXI HIGH RATE - 29 channels Total

SD PLUS FLEXI CLINIC - 29 channels Total

Number of channels: 48

No.	Channel	Maximum voltage (mV)	fc low. (Hz)	fc upp. (Hz)
1-40	EEG	$\pm 3.2, \pm 25.6, \pm 204.8, \pm 1638.4$	DC, 0.15	$0.27 * f \text{ sample}$
41-42	EMG/EP	$\pm 0.8, \pm 3.2, \pm 51.2$	0.11, 0.15, 2.2, 8.9	$0.31 * f \text{ sample}$

No.	Channel	Description
43-44	PNG	Thoracic/Abdominal Tapes
45-47	OXY	XPOD Nonin
48	MKR	$\pm 50 \mu V$ 1 Hz

EEG Amplifiers: all channels are differential gain lines; channels 1 - 24 are distinguished by the inverting input in common (G2), while the other 16 are distinguished by the possibility of selecting the inverting input connected to G2, or their corresponding inverting plug insert (black)

Noise: $0.15 \mu V_{rms}$

Input impedance: 1 GOM on amplifier

CMRR: $\geq 100 \text{ dB}$ for EEG channels and $\geq 110 \text{ dB}$ for polygraphic channels

Sampling frequency: up to 16 kHz per channel for the first 40 channels, 128 kHz for channels 41-41. Channels 43-48 are auxiliary channels for specific needs.

ELECTRICAL CHARACTERISTICS

Power supply: 5 volts DC, from Micromed BQ PCI EXPRESS interfaces or
BQ USB EXPRESS

Back electromotive force : $\pm 2.5V$

Power consumption: 600 mW

Data transfer:	digital (special differential protocol RS485)
Communication Line Isolation:	2500 volts DC, with UL1577 compliant optocouplers
Current source isolation:	3000 volts DC with DC/DC converter

MECHANICAL CHARACTERISTICS

Dimensions and weight:	19 x 13 x 2.5 cm (height); 390 g; plastic shell
Connections:	inputs with touch protection (DIN 42802); 25-pin D-Sub connector for pre-wired headphones 1 PUTP3 3-pin plastic connector for oximetry sensor block 2 bipolar connector PUTP2 for inductive strips 1 4-pin Hypertac data connector 1 RCA connector for IN/OUT triggers

ENVIRONMENTAL CONDITIONS

Operational:	temperature 10 – 40 °C Humidity 0 – 90% Pressure 70 - 105 kPa
Storage and transportation:	temperature (-20)/60°C Humidity 0 - 95% non-condensing

Pressure 50 - 105 kPa

MORPHEUS

- 16 or 24 internal channels (with 16 internal channels: 8 standard reference channels and 8 configurable as common reference or differential (EEG channels) , or monopolar channels with 24 internal channels 12 standard reference channels and 12 configurable as common reference or differential (EEG channels)), or monopolar channels (printing channels) with a bandwidth of 0.15-220 Hz
- 2 channels with direct current (for external sensors)
- 2 air pressure sensors built into DC channels, absolute and differential (only **MORPHEUS**)
- 1 serial port for oximetry sensor unit, this device does not require an external power supply); 3 channels: Plethysmogram, HR
- 3 (internal XYZ accelerometer) channels for spatial position detection (body position) (**MORPHEUS only**)
- 2 built-in circuits for XActTrace chest and abdominal inductive straps (**MORPHEUS only**)
- Variable sampling rate: sampled up to 8192 Hz per channel
- Analog-to-digital 16-bit data conversion in the 6400µV range on 24 AC channels
- Analog-to-digital 16-bit data conversion in the 1V range on two channels of direct current
- Final selection (signal-to-noise improvement) from 8192 to 512, 256, 128 Hz
- Internal microcontroller with 1024K RAM
- Graphic display for menus and viewing of received data
- Compact removable flash memory support (1 GB or more)
- internal clock
- Dual power supply: battery powered or powered by BQ USB or BQ PCI interface
- Isolation up to 2500 VAC with magnetic connectors on communication lines and DC/DC converter on power lines

Electrical cables for MORPHEUS on BQ USB or BQ PCI interface

HBC INT LTM R	length 5 m, without trigger connections, for direct communication or via TERM LTM quick disconnect system (included), without trigger connections
TERM LTM	Quick Release System for MORPHEUS Amplifiers (30cm)

Apparatus : ***MORPHEUS*** recorders

Classification : according to EN 60601-1: **ME EQUIPMENT WITH INTERNAL POWER SUPPLY, APPLIED TYPE BF** . The device must not be used in the presence of flammable anesthetics mixed with air, oxygen or nitric oxide.

Compliance : IEC 60601-1, IEC 60601-1-2, IEC 60601-2-26

Power supply : 2 x 1.5 V DC AA alkaline batteries.

Internal voltage : +3.3 3.6V

Power : 150mW

Number of channels :

- 12 (***MORPHEUS***) differential channels configurable as common reference channels: EEG channels.
- 12 (***MORPHEUS***) differential channels, bipolar configured as common reference channels or G 1: EEG/polygraphic channels
- 2 dual DC channels (for receiving data from external devices)
- oximeter at 9600 baud (no external power required), plethysmogram, SaO₂, HR
- 2 flow pressure sensor with two receiving channels (absolute and differential encoder)
- channels (built-in XYZ accelerometer) for receiving attitude data (body position sensor)
- 2 channels for inductive chest/abdominal XActTrace belts

Noise : <0.7uV r.m.s. at a sampling rate of 256 Hz

Input impedance : >10⁹ ohm

CMRR : > 100dB at 20Hz between input G 1 and all other inputs

>120dB between GND (ground) and all other inputs (IMRR)

high pass filter , 0.15Hz, 40dB/decade Hardware

low pass filter , 220Hz, 40dB/ decade 8192Hz sampling rate
per channel

Sampling rate 1024Hz, 512Hz, 256Hz, 128Hz (selectable)

ADC 16-bit for analog signals

8-bit for oximetry sensors

16-bit for internal XActTrace circuits

Dynamic range : 6400 μ V (97.6nV/value) for analog channels, 1V (15 μ V/value) for DC channels, H 2
O -sensitivity 25 μ V/cm in flow pressure channels

Impedance Control : Automatic (up to 40k Ω $\pm$ 7%): built-in generator of rectangular signals 1-16Hz.

Input Switching Built-in electronic switching of all inputs for impedance control, common
reference channel, bipolar configuration and linked to G 1.

Display and keypad Graphical, 128 x 64 pixels, 4 membrane keypad buttons and alarm

Clock Built-in, with date and alarm, can be displayed while receiving data.

Cables and connection

touch-proof connectors for connecting patient electrodes

3-pin plastic connector for daisy-chaining an oximeter

3-pin plastic connector for 2 DC channels

2 x 2-pin touch-proof connector for XActTrace (or other
connector on request).

Luer fitting for built-in air pressure transducer, metal fitting for
differential pressure transducer catheter hose.

PC connection via Micromed interface BQUSB (Direct Receive or
Download) with 4-pole plastic connector. In the process of downloading
data, the device is driven from the PC interface

- maximum speed : 4000000 bps

PC connection (Direct Data Transfer) via Bluetooth Class II radio .

- maximum speed: 460800b/s (46Kb/s)

- built-in transfer protection buffer: 16sec @ 512Hz, 32sec @ 256Hz, 64sec @ 128Hz

Communication isolation : magneto-optical connection to PC

Power isolation : built-in battery. During data download, power is supplied from the PC interface and the built-in DC / DC converter .

Data Write : Temporary to built-in RAM, semi-permanent to flash memory

Data transmission : Digital, by cable or radio channel via built-in radio module (Bluetooth )

Dimensions and weight : (2.5/3.5) H x 7.5 W x 11.4 D cm; 200g; ABS plastic body

Environmental conditions

Working : temperature 0-50°C

humidity 0-90% (excluding condensation)

Storage and transport : temperature (-20)/60° C

humidity 0-95% (excluding condensation)

SD LTM 32 EXPRESS

- **32 built-in channels for recording EEG signals (24 common reference channels, 8 configurable, i.e. used as common reference and differential channels)** with a bandwidth of 0.17-470Hz; combination of 64 devices creates an integrated system with 256 channels
- 2 additional DC channels
- 1 isolated serial port for oximeter (9600 baud rate, this device does not require an external power supply); 3 channels: plethysmogram, SaO₂, HR
- XYZ accelerometer) channels for registering the spatial position (body position) ¹
- 2 built-in inductive chest and abdomen circuits ¹
- 1 additional communication port (for communication with SD LTM STIM)
- software amplification and filtering kit
- built-in microcontroller
- graphic display for menu and windowed data acquisition
- ADC for data (pulses) 16 bits in the range of 6400 µV
- built-in memory SDHC 4GB
- built-in watch
- variable sampling rate in the range from 128Hz to 2048Hz; oversampling up to 8192Hz per channel
- two power options. Battery life over 30 hours in low power mode or power supply via interface insulation up to 2500V AC with magnetic couplings on data lines and DC / DC converter

power lines

SD LTM 64 EXPRESS

- **64 built-in channels for recording EEG signals (58 common reference channels, 6 configurable, i.e. used as common reference and differential channels) with a bandwidth of 0.17-470Hz; combination of 64 devices creates an integrated system with 256 channels**
- 2 additional DC channels
- 1 isolated serial port for oximeter (9600 baud rate, this device does not require an external power supply); 3 channels: plethysmogram, SaO₂, HR
- XYZ accelerometer) channels for registering the spatial position (body position) ¹
- 2 built-in inductive chest and abdomen circuits ¹
- 1 additional communication port (for communication with SD LTM STIM)
- software amplification and filtering kit
- built-in microcontroller
- graphic display for menu and windowed data acquisition
- ADC for data (pulses) 16 bits in the range of 6400 μ V
- built-in memory SDHC 4GB
- built-in watch
- variable sampling rate in the range from 128Hz to 2048Hz; oversampling up to 8192Hz per channel
- two power options. Battery life over 30 hours in low power mode or power supply via interface
- insulation up to 2500V AC with magnetic couplings on data lines and DC / DC converter power lines

Apparatus : SD data logger LTM 32 EXPRESS / SD LTM 64 EXPRESS
BQ USB EXPRESS / BQ PCI EXPR LTM HALF/FULL

Classification according to EN 60601-1 :

Continuous operation. The device must not be used in the presence of flammable anesthetics mixed with air, oxygen or nitric oxide.

SD LTM 32 EXPRESS : ME EQUIPMENT WITH INTERNAL POWER SUPPLY, having a means of connection to a CLASS II MAINS MAINS , TYPE BF APPLIED PART ;

SD LTM 64 EXPRESS : ME EQUIPMENT WITH INTERNAL POWER SUPPLY, having a means of connection to a CLASS II MAINS MAINS , TYPE BF APPLIED PART ;

Compliance : IEC 60601-1:2005, IEC 60601-1-2:2007, IEC 60601-2-26:2012

Power supply : SD LTM 32/64 EXPRESS : 2 x 1.5VDC AA , alkaline disposable battery, 5V DC exclusively via **BQ interface USB EXPRESS** or **BQ PCI EXPRESS**

BQ USB EXPRESS : 5V DC via compatible PC USB port according to IEC 60950 **BQ PCI EXPR LTM HALF / FULL :** 5V DC and 12V DC via PCI slot EXPRESS compatible PC according to IEC 60950

LTM HALF / FULL : 12V/5V

Input power : 2W (Max.) (BQ PCI EXPR LTM HALF / FULL) ; <1 W (typ.) / 2.5W (max.) (BQ USB EXPRESS); 750mW (SD LTM 64 EXPRESS)

Number of channels :

SD LTM 64 EXPRESS : 64 analog (16 configurable as common reference or differential)

SD LTM 32 EXPRESS : 32 analog (8 configurable as common reference or differential)

3 channels of serial port, baud rate 9600, optically isolated for connection of an oximeter

2 DC channels (for receiving data from external devices)

2 channels for thoracic/abdominal inductive belts ¹

3-axis accelerometer (body position sensor) ¹

Additional afterbirth. port 1 additional serial port for connection/exchange of data with special add. device (for example, SD LTM STIM)

Amplifiers : 32/64 channels with differential inputs distributed as: 8/16 configurable channels (differential or common reference channels) and 24/48 common reference channels

Noise : <0.5uV r . m . s . at 128Hz sample rate

Input impedance : >10 ¹² ohms at amplifier

CMRR (Common Mode Rejection Ratio) : >100dB at 50/60Hz between input G 1 and all other inputs

High Pass Prefilter : Hardware 0.15Hz, 40dB/decade **Low Pass Filter :**

Hardware 500Hz, 40dB/decade **Oversampling Rate :** 8 kHz

Sample rate : 2048, 1024, 512, 256, 128Hz

ADC : 16 bits for EEG channels, 10 bits for DC channels IN

Max. EEG input (resolution) : 6400μV (97.6nV/value)

Max. DC IN input (resolution) : 0.48/4.8V (0.46/4.69 mV/value)

Measurement accuracy : ±10%

Impedance Control : automatic (up to 40kOhm)

Display : graphic, 128 x 64 pixels

Clock : built-in, with date and alarm, can be displayed on the screen while receiving data.

Electrode box :

JB LTM 64 EXPRESS (64 channels, 6 can be used as differential or common reference, 1 common reference input, 1 GND input)

JB LTM 32 EXPRESS (32 channels: 21 head-shaped, 8 can be used as differential or common reference, 1 input - common reference, 1 input GND)

Connections : 2 x 56-pin socket box connectors.

OXI : 3-pin plastic connector for connecting an oximeter in series.

DC IN : 4-pin plastic connector for connecting DC sensors.

INTERFACE : 4-pin plastic connector for connecting to the interface (direct data reception or data download) via cable

ABDO / THOR : 2-pin keyhole connectors for recommended inductive belt models ¹

AUX : 5-pin plastic connector for connection to SD stimulator LTM STIM

standard connector USB (**BQ USB EXPRESS**) or standard PCI card EXPRESS (**BQ PCI EXPR LTM HALF / FULL**)

PC connection (direct data transfer) via radio transmission module **Bluetooth®** – class I

- maximum data transfer rate: 460800baud/sec - corresponds to 64 channels (16-bit data, sampling frequency 256Hz);

- built-in data protection buffer: 8s for 512Hz sampling, 16s for 256Hz sampling, 32s for 128Hz sampling in 64 channels

Isolation of communication channels : 2500VAC / 1min (SD LTM EXPRESS to interface) +

2500VAC / 1min (connection to PC), through magnetic couplings of data exchange lines .
2500VAC on trigger I/O lines

Power isolation : Internal battery or via PC interface and DC / DC converters (3kV DC in SD LTM EXPRESS + 3kV DC at the interface **BQ USB EXPRESS** / 1kV DC at the interface **BQ PCI EXPR LTM HALF / FULL**)

Data logging : Temporary to built-in RAM, permanently to internal micro SD - SD memory device HC 4GB.

Data transmission :	Digital, via cable or radio channel via built-in radio module
Dimensions and weight :	SD LTM 32/64 EXPRESS 4 (height) x 8 x 13 cm; 300 g; plastic case; BQ USB EXPRESS 2 (height) x 5.8 x 8.9 cm; cable length 90cm; plastic case; BQ PCI EXPR LTM HALF / FULL : Built in PC
Environmental conditions :	
Working :	Temperature 0-40°C
	Humidity 10-90% excluding condensation
	Pressure 700-1060hPa
Storage and transport :	Temperature (-20)/60°C
	Humidity 10-95% excluding condensation
	Pressure 500-1060 hPa

SD LTM 32PLUS _

- 32 internal EEG channels (24 common reference channels, 8 configurable as common reference or differential channels) with a bandwidth of 0.15-F Hz, where F depends on the sampling rate F_s ($F = 0.27 * F_s$);
- 2 auxiliary DC channels
- 1 isolated serial port for oximetry sensor box (this device does not require an external power supply); 3 channels: plethysmogram, HR
- 1 (internal XYZ accelerometer) channel for determining position in space (body position)
- 2 built-in chest and abdomen inductive belt circuit
- 1 additional communication port (for communication with SD LTM STIM)
- gain and filtering set by software
- Internal microcontroller
- Graphic display for menus and viewing received data
- A/D data conversion 24 bits in the 6400 μV range
- Support 4GB SDHC internal memory
- internal clock
- Variable sampling rate (256Hz - 16kHz) per channel
- Battery powered, run more than 24 hours in low power mode or interface powered
- Isolation up to 2500 V_{AC} with magnetic couplings on communication lines and DC / DC converter on power supply lines

SD LTM 64PLUS _

- 64 internal EEG channels (56 common reference channels, 8 configurable as common reference or differential channels) with a bandwidth of 0.15-F Hz, where F depends on the sampling rate F_s ($F = 0.27 * F_s$); combining more than 4 devices gives an integrated system

- of up to 256 channels
- 2 auxiliary DC channels
- 1 isolated serial port for oximetry sensor box (this device does not require an external power supply); 3 channels: plethysmogram, HR
- 1 (internal XYZ accelerometer) channel for determining position in space (body position)
- 2 built-in chest and abdomen inductive belt circuit
- 1 additional communication port (for communication with SD LTM STIM)
- gain and filtering set by software
- Internal microcontroller
- Graphic display for menus and viewing received data
- A/D data conversion 24 bits in the 6400 μ V range
- Support 4GB SDHC internal memory
- internal clock
- Variable sampling rate (256Hz - 16kHz) per channel
- Battery powered, run more than 24 hours in low power mode or interface powered
- Isolation up to 2500 VAC with magnetic couplings on communication lines and DC / DC converter on power supply lines

Classification according to EN 60601-1:

The device must not be used in the presence of flammable anesthetics mixed with air, oxygen or nitric oxide.

SD LTM 32 PLUS: ME EQUIPMENT with an INTERNAL POWER SUPPLY having means of connection to a CLASS II SUPPLY MAINS , TYPE BF APPLIED PART;

SD LTM 64 PLUS: ME EQUIPMENT with an INTERNAL POWER SUPPLY having means of connection to a CLASS II SUPPLY MAINS , TYPE BF APPLIED PART;

Compliance with standards: IEC 60601-1, IEC 60601-1-2, IEC 60601-2-26

Power supply:

SD LTM 32 PLUS: 2 x 1.5 VDC AA, non-rechargeable alkaline battery and 5 VDC exclusively from BQ USB PLUS, BQ USB MULTI or BOX interface TERMINAL MULTI . Rechargeable internal lithium battery.

SD LTM 64 PLUS: 2 x 1.5 VDC AA, non-rechargeable alkaline battery and 5 VDC exclusively from BQ USB PLUS, BQ USB MULTI or BOX interface TERMINAL MULTI . Rechargeable internal lithium battery.

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Input power: 1.5 W (SD LTM 64 PLUS);

Number of channels:

SD LTM 32 PLUS: 32 analog (8 configurable as common reference or differential)

SD LTM 64 PLUS : 64 analog (8 configurable as common reference or differential).

3 channels per serial port, 9600 baud, optically isolated for oximetry sensor block

2 DC channels (for collecting data from external devices)

2 channels for chest / abdomen induction tape

Triaxial accelerometer (body position sensor)

104

Auxiliary serial port 1 auxiliary serial port for connection/communication with certain auxiliary devices (e.g. SD LTM STIM)

Amplifiers:

32 channels with differential inputs divided by 8, configurable as differential or common reference and 24 with fixed common

64 channels with differential inputs divided by 8, configurable as differential or common reference and 56 with fixed common

Reference

Noise: 0.2 μ Vrms @ 256Hz sample rate

Input impedance: >109 ohm (amplifier)

CMRR: >100 dB @ 50 Hz between G1 and all other inputs

High Pass Filter: DC, 0.15 Hz

Low pass filter: $0.27 * f_{\text{sampling}}$

Sample rate: 16384, 8192, 4096, 2048, 1024, 512, 256Hz

Analog to digital conversion: 24 bits for EEG channels, 10 bits for DC IN channels

Maximum EEG input (resolution): full scale range: $\pm 3.2, \pm 25.6, \pm 204.8, \pm 1638.4$ mV resolution up to 44 nV (24 bit ADC)

Max DC Input (Resolution): 4.8V (4.69mV/digit)

Measurement accuracy: $\pm 10\%$

Impedance Check: Automatic

Display: graphic, 128x64 pixels

Clock: internal, permanent with date and alarm

Jackbox JB LTM 64 PLUS (64 channels, 6 can be used as differential or common reference, 1 common reference input, 1 GND input)

JB LTM 32 PLUS (32 channels: 21 head, 8 can be used as differential or common reference, 1 common reference input, 1 GND input)

Connections: 2x56pin Connectors for jackbox.

OXI: 3-pin plastic connector for daisy-chaining to the oximetry sensor box

DC IN: 4-pin plastic connector for connecting to DC sensors.

INTERFACE: 4-pin plastic connector for connecting to the interface (straight data collection or download) via cable

ABDO / THOR: Keyed 2-pin connectors for approved inductive strips

AUX: 5-pin plastic connector for connection to SD LTM STIM stimulator

Standard USB connector (for BQ USB PLUS and BQ USB MULTI)

Communication with PC via wired network (TCP/IP)

Communication with PC via Wi-Fi radio transmission module (IEEE 802.11 b/g Wireless LAN) with secure encryption and roaming capability: configurable wireless channels

Communication line isolation: 2500 VAC / 1 minute (SD LTM PLUS to interface) + 2500 VAC / 1 minute (PC interface) via magnetic couplings on communication lines.

2500 VAC on trigger I/O lines

Power supply isolation: internal battery or from PC interface via DC/DC

converters (3 kV DC to SD LTM PLUS + 3 kV DC to BQ USB PLUS interface + 5 kV DC to BQ USB MULTI / BOX TERMINAL MULTI)

interface)

Data recording: temporarily in internal memory. Constantly on the internal Micro SDSD HC supports 8GB Scandisk Industrial.

Data transmission: digital, cable or radio via built-in Wi-Fi module

Dimensions and weight: SD LTM 32/64 PLUS 4 (height) x 8 x 13 cm; 300 g; plastic case;

Environmental conditions:

Working temperature 0-40°C

Humidity 10-90% non-condensing

Pressure 700-1060 hPa

Storage and transport: Temperature (-20) / 60 °C.

Humidity 10-95% non-condensing

Pressure 500-1060 hPa

Isolating transformer TR 50

Classification according to Directive 93/42 / CEE: class I device

Classification according to EN 60601-1

Electrical safety : **ME EQUIPMENT CLASS I, APPLIED TYPE B**

The device is not protected against ingress of water or dust

The device is suitable with continuous operation mode

The device is not suitable for use in the presence of flammable anesthetic mixtures with air, oxygen or nitrogen.

Storage conditions: temperature -20 / 60 ° C, humidity 0-95% non-condensing

Operating conditions: temperature 0-45 ° C, humidity 20-95% non-condensing

Power supply: 210 ~ 240VAC, 50Hz

Maximum power consumption: 500 VA

Mains connection: VDE male type (IEC320/C14)

(external power cable, minimum power 10A, maximum length 180cm.

VDE socket (IEC320 / C13)

Fuse: 2x 4A T

Outputs: 6 SCHUKO connectors (CEE 7/7)

220~250V, 50Hz

parallel connection (maximum output power 500 W; output power is divided between outputs by load)

Thermal protection: 110°C

Conforms to IEC 60989 and EN 60601-1 safety standard

Isolation transformer TOR 500

Classification according to Directive 93/42 / CEE: class I device

Classification according to EN 60601-1

Electrical safety : ME EQUIPMENT CLASS I, APPLIED TYPE B

The device is not protected against ingress of water or dust

The device is suitable for continuous operation

The device is not suitable for use in the presence of flammable anesthetic mixtures with air, oxygen or nitrogen.

oxide.

Storage conditions: temperature -20 / 60 ° C, humidity 0-95% non-condensing

Operating conditions: temperature 0-45 ° C, humidity 20-95% non-condensing

Power supply: 210 ~ 240VAC, 50Hz

Maximum power consumption: 500 VA

Mains connection: VDE male type (IEC320/C14)

(external power cable, minimum power 10A, maximum length 180cm.

VDE socket (IEC320 / C13) for connection according to the sockets used by the user
a country)

Fuse: 2x 4A T (delay)

Outputs: 6 SCHUKO connectors (CEE 7/7)

220~250V, 50Hz

parallel connection (maximum output power 500W; output power is divided between
load outputs)

Thermal protection: 110°C

Conforms to IEC 60989 and EN 60601-1 safety standard

LTM TP BUTTON

It is an accessory for devices of the SD LTM EXPRESS family and must be used for the intended purpose of the devices to which it is connected.

- Description: The LTM BUTTON TP is the patient button for the heads of the SD LTM EXPRESS family. Allows the user to insert a marker (dot on the graph) on an EEG trace received by the device it is connected to.
- Main material: self-extinguishing PVC
- Dimensions: cable length 30cm
- Reusable: yes, non-sterile and non-sterilizable device
- latex-free: yes

Applicable regulation

Read / Guide

Directive 93/42/EEC

Legislative Decree 46/1997 and subsequent amendments

Technical standards

EN 60601-1

classification

second set / product type / grade notes

Directive 93/42 / EEC Rule II CERTIFICATION PROCEDURE:

Appendix VII: Manufacturer's Declaration of Conformity

PRODUCT CODE 380.040.0012

TYPE OF DEVICE

(classification according to GMDN nomenclature) N.D.

CND CODE Z 12101003 - EEG HOLTER RECORDERS

REPORT CODE 24182

Intended use

MORPH BUTTON is an accessory for MORPHEUS devices and must be used for the intended purpose of the devices to which it is connected.

Product Description

MORPH BUTTON - Allows the user to insert a marker (dot on the graph) on the EEG trace received by the device it is connected to.

- Main material: ABS plastic
- Measurements: button 3.1 x 5.6 x 3 cm; weight: 100 g; cable length 1m.
- Reuse: reusable, non-sterile and non-sterilizable device. Cleaning with a damp cloth (water and mild soap solution); disinfection with 1:2 diluted isopropyl alcohol.
- latex-free: yes

Applicable regulation

Directive 93/42/ EEC and subsequent amendments transposed into Legislative Decree 46/1997 and subsequent amendments

Technical standards

CEI EN 60601-1

classification

second class/type notes

Directive 93/42/ EEC I CERTIFICATION PROCEDURE:

Annex VII : manufacturer's declaration of conformity

Cortical Stimulator SD LTM STIM

Electrical specification:

Stimulus waveform: positive or negative impulse. Positive and negative pulse with selectable first edge and separate independent duration. Intensity: 0.1 to 15 mA adjustable in 0.1 mA steps

Single pulse parameters:

Positive Pulse Width: 50 to 3000 μ s (in 50 μ s increments)

Negative Pulse Width: 50 to 3000 μ s (in 50 μ s increments)

Positive-negative distance: 25 - 100 μ s

Continuous stimulation (frequency <20 Hz):

Frequency: 0.1 to 19.9 Hz

Sequence duration limit: 1 to 500 s.

Continuous stimulation (frequency > 20 Hz):

Frequency: 20 to 300 Hz

Sequence duration limit: 1 to 15 s.

Continuous stimulation of the line (frequency <20 Hz):

Sequence duration limit: 1 to 15 s.

Line frequency: 0.1 to 19.9 Hz

Pulse period: 3.3 to 10000 ms

Pulse frequency: 0.1 to 300 Hz

Line Distance: ≥ 10 ms

Maximum stimulation output voltage: 26V

The higher impedance at which 15 mA stimuli can be delivered is about 1.7 k Ω , at higher impedances the current decreases. The error in the value of stimulus intensity does not exceed 10%.

for patient impedances below 1.7 k Ω and for currents greater than 0.5 mA (for current settings for which pacing delivery is possible).

Stimulating electrode selector

Independent 64-pole selector for 64-channel electrodes

Isolation of stimulating electrodes via Photo MOS and reed relays

Power supply

Via a cable connected to the SD LTM 64 EXPRESS recorder and 3 AAA 1.5 V batteries located inside the stimulator. Alkaline batteries are recommended for high power electronic devices (such as cameras). Both power supplies must be present. The yellow STIMULUS LED is used to alert the user to the state of charge of the battery.

Classification: according to EN 60601-1: **ME EQUIPMENT WITH INTERNAL POWER SUPPLY, APPLIED PART TYPE BF**. To meet the standard for this type of classification, the stimulator must only be used with the SD LTM 64 EXPRESS.

Power supply: SD LTM 64 EXPRESS and 3 AAA 1.5V batteries (both required). The SD LTM STIM stimulator and the SD LTM 64 EXPRESS amplifier together form a medical electrical system.

Insulation: to PC over 3000V.

IP classification: IP20 (Protection against penetration of solid external objects with a thickness of 12.5mm.

Ø and above; not protected against penetration of liquids)

Size: 9 x 10 x 4 cm.

Environmental conditions

Operating: temperature 0-40°C, humidity 10-90% non-condensing

Storage/transport: Temperature (-20) -60°C, humidity 10-95% non-condensing.

Output specifications: Electrical stimulation:

Electrical stimulation with constant current.

Controlled safety limits (charge per phase less than 20 μC)

Stimulus Intensity: 0 to 15 mA in 0.1 mA steps

Stimulus duration: from 50 μs to 3000 μs in 1 μs steps

Stimulus repetition: single impulse; 0.1 Hz to 300 Hz in 0.1 Hz steps

Stimulus mode:

Monophasic (positive or negative)

Two-phase (with independently adjustable main and reverse phases)

Maximum stimulus intensity: 15 mA (load up to 1.7 k Ω)

Stimulator max. operating voltage: 26 IN

Maximum charge per phase: 20 μC

Block for connecting oximetric sensors (the system does not measure, but uses sensors to monitor the condition).

Refresh rate value: 1 Hz

Accuracy (± 1 standard deviation)

70-100% ± 2 digits with Finger Clip Sensor

accuracy not specified for values <70%

heart rate:

range 18-300 bpm

accuracy $\pm 3\% \pm 1$ digit

update rate: 1 Hz

Plethysmogram:

graphic representation

display the update rate associated with the timebase

Sampling frequency:

75 Hz at origin

64 samples stored per second

Photostimulator LED with bracket

Classification according to IEC 60601-1: **CLASS I ME EQUIPMENT, TYPE BF APPLIED PART**

Classification according to IEC 62471: Risk group I

Storage and transportation conditions: temperature -20 - 60 ° C

humidity 10 - 95% non-condensing

pressure 500 hPa - 1060 hPa

Operating conditions: temperature 0 - 30 ° C

humidity 20 - 90% non-condensing

pressure 700 - 1050 hPa

Input voltage: 100-240VAC, 50/60Hz

Input power: max. 15 VA (at 60 Hz)

LED photo equipment

Lamp supply voltage: DC 12V

Flash period: 2ms

Flash frequency: 0... 60 Hz

Light discharge hole: 123 x 31mm

Light Intensity:

Measuring distance 10 cm: 30 Hz: 1.0 klux

60 Hz: 1.9 klux

Measuring distance 20 cm: 30 Hz: 0.6 klux

60 Hz: 1.1 klux

LED used: ML5FW13H-CEF

LED brightness: 4000 mcd

LED Beam Angle: 20°

Dimensions: 173 × 55 × 70 mm

LED PHOTO POWER

Fuses: 2 × T 1AH, 250 V, 5 × 20 mm

Trigger input specifications: 0-5V positive

Threshold trigger input (on): ≥ 2.7 V and ≥ 1.3 mA

Threshold trigger input (off): ≤ 1.2 V and ≤ 0.3 mA

Dimensions: 153.5 × 50 × 89.5mm

LED PHOTIC System User Manual 12

Power cable (for EU market): Schuko cable (CEE 7/7) - UEC320 C13

Voltage: 250V

Current: 10A

Conductor cross section: 0.75 mm²

Power cord (for US market): NEMA Cable 5-15P - IEC320 C13

Voltage: 125V

Current: 15A

Conductor section: 16AWG

*Note: Technical specification for optional parts and accessories is given in the Appendices to this document.

1.2.5 The design of the EEG system should ensure the optimal use of standard and reusable design solutions, a rationally limited range of products, grades and range of materials.

1.2.6 The design should provide maximum ease of maintenance. All surfaces of the EEG system must be accessible and suitable for cleaning and washing with detergents.

1.2.7 The design of the components of the EEG system should prevent the possibility of their incorrect assembly (connection). The mechanical parts to be joined must be tightly connected to each other without gaps.

1.2.8 The degree of protection of electrical equipment is indicated in the table of technical characteristics

1.2.9 For modes of non-automatic analysis of EEG signals in the EEG system, amplitude and time stamp calibrators must be provided. The relative error of registration of the calibration signal should be no more than $\pm 5\%$.

The calibration signal in this case must be entered in front of all elements of the alternating current coupling, which determine the time constant, and in front of all regulators accessible to the operator (with the exception of switching devices of leads).

1.2.10 Operating conditions of the EEG system

1.2.11.1 The EEG system is suitable for operation under the following conditions:

Permissible ambient temperature during operation: from plus 10 to plus 35 °C, humidity: from 25 to 95% at plus 25 °C, atmospheric pressure: from 86 to 106.7 kPa .

1.2.11.2 The packaged EEG system must retain its appearance and operability after exposure to multiple impact loads with a peak impact acceleration of not more than 9 g with an impact acceleration duration of 10 ... 15 μ s, or transportation along dirt roads (at least 200 km) .

1.2.11.3 Cables tensioned with a force of 1 N must withstand at least 500 bends of $\pm 90^\circ$ at the point of their exit from the housing.

1.2.12 The EEG system should be designed so that its electronic part is maximally adapted to the use of automatic assembly.

1.2.13 The EEG system must contain a minimum number of electronic part settings, but not to the detriment of its technical parameters.

1.2.14 Equipment from the EEG system must have a smooth, uniform and one-color outer surface of the body, without burrs, cracks, stains, dirt, scratches, corrosion, pits, dents and chips, corresponding to the approved standard samples.

The roughness of the surfaces accessible to touch by the operator should be no higher than $R_{\text{aut}} 6.3$ microns

1.2.15 Controls must be marked, signs and symbols must be clear.

1.2.16 The surface of the parts intended for applying protective coatings must be clean, even and smooth: oil stains or other contamination, waviness, large scratches and rough roughness are not allowed.

The EEG system must be resistant to short circuits and open circuits .

1.2.17 Reliability indicators

1.2.17.1 The established failure-free operating time shall be at least 1,500 hours.

The mean time between failures must be at least 5,000 hours.

Failure of the EEG system is a violation of its operable state associated with the failure of any component, resulting in a deviation of the parameters beyond the limits established in these technical specifications, if, in order to restore the operable state, it is necessary to replace or repair any component.

1.2.17.2 The average recovery time of a working state when replacing a defective part (element) with a serviceable one should be no more than 2 hours.

1.2.17.3 The total specified service life of the EEG system is 5 years.

1.2.17.4 EEG system technical utilization coefficient – not less than 0.90.

1.2.17.4 In terms of the reliability of the transmitted information, the EEG system must comply with :

- the probability of team transformation should not exceed 10^{-14} ;
- probability of forming a false team: 10^{-12} ;
- probability of refusal to execute the sent command (it is allowed to repeat the transmission up to 5 times): no more than 10^{-10} ;
- probability of information loss during sporadic transmission (repetition of transmission up to 5 times is allowed): no more than 10^{-8} .

1.2.17.5 The probability of a failure leading to a false operation of the modular system for 1,000 hours of operation is not more than 0.01.

The coefficient of undetected errors is not higher than 10^{-4} .

1.2.18 The service life of the electrodes is 1 year.

1.2.19 The EEG system software must have:

- noise immunity, i.e., it eliminates the possibility of looping and falling into deadlock situations, the ability to function correctly in case of failures, failures of a part of computing facilities, personnel errors, and also provides automatic restart when power is restored after a shutdown without issuing false signals and control actions;

- performance, i.e., minimum restart time, reactions to external events, minimum loss of time for data protection;

- adaptability and development prospects, i.e., the ease of adapting programs to changes or extensions of user tasks without compromising other indicators.

1.2.20 The duration of daily maintenance should be no more than 1 hour, monthly - no more than 5 hours, annual - no more than 10 hours. The following information is for daily maintenance at the location of the system, monthly and annual maintenance should be contacted by the manufacturer or his authorized persons.

ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL STUDIES MICROMED , PERFORMED BY: SD LTM 32 EXPRESS

ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL STUDIES MICROMED , IN VERSIONS: SD LTM 64 EXPRESS

SD LTM 32 EXPRESS and SD LTM 64 EXPRESS include two power supply systems that increase flexibility in working with it. In amplifier mode « amplifier mode ».

be used for a short period of time to allow the patient to move around the ward. In "ambulatory mode" the batteries are used as the primary power source since the cable is not connected.

The power supply at the interface has safe characteristics. Interface output isolated from PC USB / PCI connector EXPRESS via converter DC / DC , and the trigger input / output is isolated from the USB / PCI connector at the same time EXPRESS and communication lines with an amplifier. To **SD amplifier LTM 32 EXPRESS** a 5V power line is supplied, which is converted at the next stage of isolation.

DO NOT USE A POWER DIFFERENT FROM THE INSTRUCTIONS (BQ USB EXPRESS / BQ PCI EXPR LTM HALF / FULL AND 1.5V AA BATTERIES), FOR THE OPERATION OF THE AMPLIFIER: OTHERWISE SERIOUS DAMAGE TO THE INSTRUMENT MAY OCCUR WITH A RISK TO THE OPERATOR AND THE PATIENT.

Do not connect to the interface via cable, **SD LTM 32 EXPRESS and SD LTM 64 EXPRESS** use two 1.5V batteries for their operation. Two new batteries can work for 24 hours of continuous recording at 128Hz sampling rate (Duracell batteries recommended) Ultra 1.5V AA); in standby mode, they can work for several months.

During data reception, the following symbols appear on the display:

 the batteries are good (the scale decreases as the level of charge of the batteries decreases);

 The batteries are dead and must be replaced.

If the batteries are low and the power cable is not connected, the recorder emits a sequence of 3 short beeps to warn of this and the red LED indicator flashes constantly. Batteries must be replaced.

If the batteries are replaced while receiving data, and the device is turned on, the user is prompted to decide whether to continue with the previous data receive or start receiving data again.

If the device is used as an ambulatory recorder without a cable connected to a PC, the batteries should be replaced as soon as possible if the battery status indicators (LED, display indicator, PC screen messages) come on and alert the user. If the recorder operates on incompletely charged batteries and data is not received for a period of time according to the device settings, then data reception will stop before the examination is completed, or an erroneous recording will occur (for example, memory sections are not available to record a new pattern). See the SAFETY PRECAUTIONS section and follow the operations in the instruction manual.

IMPORTANT: always check the condition of the batteries before turning on the recorder.
Do not use batteries that are not fully charged.

Remove the batteries if the device is not expected to be used for an extended period of time.

AA alkaline batteries to ensure the performance of the recorder. Follow the manufacturer's instructions for using and disposing of batteries.

Pay attention to the polarity of the batteries! Before inserting the batteries, check the correct position marked on the base of the battery holder.

To prevent damage to the device when replacing batteries, follow the following procedure for installing them.

Correct: install the battery with correct polarity by inserting the right side of the battery into the holder first, then press the other side of the battery with your finger so that the left side of the battery takes its place in the holder.

Incorrect : The battery is installed in such a way that the left side of the battery is inserted first, and then the right side

Recording duration

Pattern files are the same size and fill the entire amount of memory except for sectors that are reserved for re-programming the device. The file size can only be changed by new segmentation/partitioning of the memory card in PARTITION mode IN USE in the main menu in the program.

The received EEG data is written to a file, which is determined by the partitioning of the memory into sections: the maximum length of the record is related to the number of available sections, the memory size and the sampling frequency (sampling rate), the number of channels and the type of

Please note that the registration time marked with * exceeds the performance of the batteries. Therefore, batteries need to be replaced while data is being received. It is possible to save the received data in the same file after replacing the batteries if the pause (break) duration is less than 5 minutes. The time elapsed from the moment data reception was stopped until the moment it was restarted is filled with isoelectric (zero signal). Using a wireless connection increases power consumption and battery life may be shortened.

Memory size	Mode	Channels and sample rate	Duration
4 GB	1 partition, uncompressed data	32 channel 128Hz	Approximately 149 h *
		32 channels 128Hz + 2 channel DC AUX 128Hz + oximeter data (SaO ₂ , HR , plethysmogram)	Approximately 143 h *
		32 channel 256Hz	Approximately 74.5 h *
		64 channel 128Hz	Approximately 74.5 h *
		64 channels 128Hz+	Approximately 73 h *
		2 channel DC AUX 128Hz +	
		64 channel 256Hz	Approximately 39 h *
		64 channel 512Hz	Approximately 19.5 h
		64 channel 1024Hz	Approximately 9.5 h
		64 channel 2048Hz	Approximately 4.5 h
		64 channel 2048Hz+	Approximately 4.5 hours
		2 channel DC AUX 128Hz + oximeter data (SaO ₂ , HR , plethysmogram)	

Standard

$duration [h] \approx ((size [MB] \div 1048576 \div 500 \times 1024) / n \div partition) \div 512 \times 64 ((n \div channels_EEG \times samp.\ freq._EEG \div 128^{(a)}) \div 2 \div 66^{(b)}) \div 3600$

samp . freq . _ EEG - standard sample rate per channel

n ° channels _ EEG number of channels that are used by the data receiving program (from 1 to 64)

(a) 128 must be added if DC channel data reception is enabled AUX ;

(b) 66 should be added if oximeter and plethysmogram data are obtained.

Device connectors

Assignment of pins in the connector from the side of the cable connection.

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

Maintenance and updating of hardware and software

The microcontroller and CPLD memory in the recorder can be reprogrammed directly by the user in the SystemPlus program EVOLUTION .

the SD firmware **LTM 32 EXPRESS** and **SD LTM 64 EXPRESS** through the cable connected to the interface. The file used for the update (NEWLTMB1 . DAT) is contained in the SystemPlus installer EVOLUTION .

The progress bar shows the status of the copy operation; at the end of copying, the process of updating the program begins. The on-screen counter shows the number of remaining blocks that need to be reprogrammed. Then a corresponding message appears informing the user about the results of the operation. If everything is in order, the device will automatically restart with the new firmware. If not, the recorder must be reprogrammed.

Automatic testing

The recorder firmware contains several applications for testing the health of the analog and digital blocks of the device.

Analog block tests are available in the impedance control menu by pressing the **UP button** . The results are output in the following order; to move from one test to another test, press the **UP** arrow button , to return to the previous test, press **DW** , to exit test mode, press the **OK button** (**ESC**).

- Gain inverting and non-inverting inputs using a test signal (16Hz, 5mV square wave).
- Alignment of inverting and non-inverting inputs (compared to channel 3) using a test signal (16Hz, 5mV square wave).
- ☐ Biasing channels with inputs connected to ground using 50kΩ resistors.
- Background noise at inputs connected to ground through 50kΩ resistors.
- ☐ External signal amplification (it is required to apply an external signal with a voltage of 3200 μV to the inputs of the recorder).

All results are displayed in columns whose height depends on the value measured in each channel. The Y- axis is the scale for the current test values (the upper right corner of the display shows the full size value). Maximum and minimum values on the scale The y -axes represent the limit values: the recorder works properly if all tests give results that fall within the range of standard values.

The digital block test can be performed at the time of startup by simultaneously pressing the **ON buttons** And **OK** . The test function tests all memory by monitoring data and addresses. At the end of the test, the test result is displayed: if the test fails, the message also contains the section of memory where the first error was found.

Tests using hotel tools

MM 10 KRES cable set or a connector with 10kΩ resistors, which are connected to the input of the

signals at the input to the recorder.


A voltage divider ratio of 1:1000 should be used.

System Plus Software **EVOLUTION** or **Device** utility **test** should be used during the measurement process, so it is more convenient to use the device in the “electrode box” mode.

Device utility **test** follow the instructions and options that the program requires.

The setting specified in the following tests may differ from the settings required in the **Device Utility test**.

Test CMRR– EEG -channels

Rice. 10 connection diagram for test execution CMRR

The common mode rejection ratio (CMRR) test is performed by setting a balanced sine wave at 20Hz, 2V amplitude, zero offset for each pair of inverting and non-inverting inputs. Set “ connect ” (data receiving program) for 64 common reference channels (G 2) and set the sampling rate to 256Hz. Turn on the data reception mode.

View data in System Plus by setting up editing for 64 common standard channels in open mode (no software filters) and no notch filter.

Based on this value, you can calculate the CMRR in dB using the formula

$$CMRR = 20 \log_{10} \left(\frac{2}{v} \right)$$

where v is the doubled signal amplitude of each channel.

This parameter must be greater than or equal to 100dB.

CMRR coefficient value and the measured amplitude is shown below.

μV PP	CMRR (VIN
1	126
2	120
3	116.5
4	114
5	112
6	110.5
7	109.1

mkV PP	CMRR (VIN=2V)
8	108
9	106,9
10	106
11	105,2
12	104,4
13	103,7
14	103,1

mkV PP	CMRR (VIN=2V)
15	102,5
16	101,9
17	101,4
18	100,9
19	100,4
20	100
21	99,6

Background noise test

Shorting inverting and non-inverting inputs to ground.

Set the sampling rate to 128Hz and the printing channels to differential mode. Perform data logging.

View data in System Plus , set up editing to view all recording channels in clear mode (no software filters) and no notch filter.

Check that the background noise is below 5μV PP for 0.5s.

Set the signal generator to output a 100mV pp square wave at 10Hz. Use a balanced voltage divider input, apply a signal to the recorder with 1000 times attenuation (10Hz input becomes 100μV pp). Connect the positive pulse output of the voltage divider to the non-inverting input (channels 1 to 64) and the negative pulse output of the voltage divider to the inverting input (channels 25...30 and G 2) and ground to GND .

Set up the program in System Plus to view all printing channels in differential mode. Set the filter to the widest frequency band. Turn on data reception and check that the signal amplitude is $100 \pm 10 \mu V$.

AUX DC Channel Test

pp square wave at 8Hz. Set channels DC 1/ DC 2 to 4.8V.


Apply a signal to DC 1 and DC 2 with a unity gain divider. Receive data.

Check that the received signal is $2 \pm 0.2V$

Impedance control

Connect 10kΩ resistors between the inverting and non-inverting input and ground for all 64 channels. Use an impedance checker to make sure the reading is correct.

Input switch test

RICE. 12 Wiring diagram for switch test

pp square wave at 10Hz. Use a balanced voltage divider, apply a 1000 times attenuated signal to the input of the recorder (input signal must be 100mV pp , frequency 10Hz). Connect the positive pulse output of the divider to the non-inverting input of channels 25 to 30 and to the non-inverting input of all channels. Connect the center point of the divider to the inverting input of the common reference channels (G 2), and the negative pulse output to GND (G 1).

In System Plus EVOLUTION set up three connections that configure the channels in differential mode, common reference mode, or linked to G 1. Select the widest bandwidth.

Enable data reception using each connection (in case of data reception by a program associated with the recorder). Apply a square wave pulse to the positive and negative pulse inputs of the configurable channels from the positive output

of the divider).

Check that the amplitude of the signals in the configurable channels is:

- 0V if the channel is configured as differential;
- 50 □ 10μV if the channel is configured as a common reference (the input of negative pulses is connected to G 2);
- 100 □ 10μV if the channel is tied to G 1.

ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL STUDIES MICROMED, PERFORMED BY: SD LTM 32 PLUS AND ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL RESEARCH MICROMED , VERSION: SD LTM 64 PLUS :

Specifications of **SD Plus Amplifier** can be monitored with special tests using laboratory instruments and an admission program. A package called **MICROMED DEVICE TEST** is also available , containing a special SD PLUS TEST program that automates the control procedure and test report generation.

For the controls described below, it is recommended to equip the connectors with 10 kΩ plug resistors at the amplifier inputs. Moreover, it is necessary to use a voltage divider function generator to generate proper signals to the amplifier.

It is necessary to have a function generator with an external trigger that allows you to adjust the number of signal cycles.

In all tests below, the external trigger of the signal generator must be connected to the red RCA connector located on the communication interface between the PC and the amplifier or on the hollow punch of the interface board. It is necessary to set the number of cycles equal to five with a signal produced by the generator.

A voltage divider with a ratio of 1:1000 must be used.

Measurements must be made using the SystemPlus EVOLUTION software.


Fig. 12. Wiring Diagram for CMRR Test (Common Mode Rejection Ratio)

The CMRR test is performed by sending a balanced sinusoidal signal (make sure the device offset is zero), at a frequency of 50 Hz and an amplitude of 1.6 V, to all EEG channels, including G2 (negative input pulse of the common reference channels) and G1 (control back electromotive force point). Turn on the EEG acquisition program mode ALL CHANNELS and set the sampling rate to 512 Hz.

Set the filters to the maximum bandwidth of high frequencies 0 Hz, low frequencies 1000 Hz, capture window 50 $\mu\text{V}/\text{cm}$.

Turn off the notch filter and start receiving.

This information allows us to calculate

$$\text{CMRR} = 20 \log_{10} \left(\frac{1.6}{V} \right)$$

where V is the signal span of each channel

This parameter must be greater than or equal to 125 dB.

The correspondence between some CMRR values and the measured amplitude is shown below.

μVPP	CMRR (@ $V_{\text{in}}=1.6\text{V}$)	μVPP	CMRR (@ $V_{\text{in}}=1.6\text{V}$)	μVPP	CMRR (@ $V_{\text{in}}=1.6\text{V}$)
0.1	144	2	118	12	102
0.2	138	3	114	13	101
0.3	134.5	4	112	14	101
0.4	132	5	110	15	100
0.5	130	6	108	16	100
0.6	128.5	7	107	17	99
0.7	127	8	106	18	98
0.8	126	9	104	19	98
0.9	125	10	104	20	98
1	124	11	103	21	97

Table 6 Correspondence between CMRR and measured voltage.

CMRR test – polygraphic channels

The CMRR check is carried out by sending a balanced sinusoidal signal, with a frequency of 50 Hz and an amplitude of 1.6 V, to all 11 printing channels. Set the EEG acquisition program mode to C21P and set the frequency of the 11 polygraphic channels to 1024 Hz.

Set the filters to the maximum bandwidth.

Output the amplitude of the signal from the video.

This information allows you to calculate the CMRR per dB using the formula and values in Table 6: this parameter must be greater than or equal to 125 dB.

Background noise test – EEG channels

Short circuit all 21 EEG inputs and G2 to G1.


Set the sampling rate to 256 Hz and the filters to a bandwidth of 0.3 - 30 Hz. Start a receive and verify that the background noise is below 0.9 μV peak-to-peak (0.15 μV rms) within a 0.5 s time window.

Background noise test – polygraphic channels

Short all inputs (positive and negative) of 11 printing channels to G1.

Set the sampling rate to 256 Hz and the filters to a bandwidth of 0.3 - 30 Hz. Start a receive and verify that the background noise is below 0.9 μV peak-to-peak (0.15 μV rms) within a 0.5 s time window.

Proper amplification test - EEG channels

Pic. 11. Amplifier Test Connection Diagram

Set the signal generator to deliver a 5 Vp-p, 10 Hz square wave. Use a balanced voltage divider to send the generated signal attenuated by a factor of 1000 (yielding 5000 μ V peak-to-peak at 1 Hz) to the amplifier's printing inputs. Connect the positive output of the divider to the positive input of the amplifiers (channels 1-24), the negative output of the divider to G2, and ground to G1.

In the BQ software, select the program showing all EEG channels with reference to G2. In this program, you need to set the filters to the widest bandwidth: 1024 Hz. Turn on the "all channels" reception mode and select the 6400 μ V reception window. Verify that the observed amplitude is 5000 $\pm$ 500 μ V.

Background noise test – polygraphic channels

Set the signal generator to deliver a 5 Vp-p, 10 Hz square wave. Use a balanced voltage divider to send the generated signal attenuated by a factor of 1000 (yielding 5000 μ V peak-to-peak at 10 Hz) to the amplifier's printing inputs. Connect the positive output of the divider to the positive input of the amplifiers (channels 25 - 40), the negative output of the divider to the negative output of the amplifier, and ground to G1.

In the receiving software, select the program that shows all printing channels. In this program, you need to set the filters to the widest bandwidth. Set the sampling rate to 128 Hz. Turn on the "all channels" reception mode and select the 6400 μ V reception window. Verify that the observed amplitude is 5000 $\pm$ 500 μ V.

High bandwidth Vp channel test

Set the function generator to deliver a 5 Vp-p sine wave at 10 Hz. Use a balanced voltage divider to send 1000 times the generated signal to the EP1 and EP2 inputs. Verify that the observed amplitude is 5000 $\pm$ 500 μ V.

Impedance check

Connect 10 k Ω resistors between the inputs of 40 EEG/EP channels and ground (G1). Verify that the readings are correct using an impedance checker. Carry out control of inputs EP1 and EP2 in the same way.

Autodiagnostic program

The SD Plus Family Amplifier Self-Diagnostic Program is available upon request .

The program allows the user to test background noise and correct gain, measure input impedance, and function check the built-in switches to change filters and receive windows. In addition, matching amplifier input connectors can be supplied to speed up tests. The program requires the use of a function generator.

ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL STUDIES MICROMED , PERFORMED BY: SD PLUS FLEXI CLINIC , ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL STUDIES MICROMED , PERFORMED BY: SD PLUS FLEXI , ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL STUDIES MICROMED , VERSION: SD PLUS HIGH RATE, ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL STUDIES MICROMED , PERFORMED BY: SD PLUS RESEARCH:

Check the box, accessories, and recording cables for accidental shipping damage. Possible damage must be examined and reported to your local Micromed representative before using the instrument.

Before cleaning the system, make sure that the power supply is turned off. Do not use any abrasive materials or solvents for cleaning.

Clean the device with a damp (not wet) cloth on the outer sides of the case. Cables should also be cleaned with a damp cloth. For hygiene reasons, it is necessary to replace the electrodes after each use.

If the device has been dropped from a height not exceeding 1 m, before starting work on the patient, it is necessary to verify the normal functioning of the device using an impedance test check (CHECK) and calibration resistors. If dropped from heights above this value, the device must be checked by Micromed personnel .

Do not open the system cover or expose it to abnormal operating conditions (high temperature and excessive humidity).

Avoid getting any materials, especially liquids, on the system components, in which case, do not use the system under any circumstances, and contact our company representative in your area for assistance.

Do not repair or replace any part of the system yourself. In the event of a malfunction, contact qualified personnel. In case of potentially dangerous situations, it is necessary to check the head module, even if it seems to be operational.

The microcontroller and programmable logic in the amplifier can be reprogrammed by the user directly through the SystemPlus program EVOLUTION .

The device settings screen (see *SETUP PROCESS* on page 20) has a button that allows you to flash the SD firmware Plus via connection cable and interface. The file used for the update (SDPlus . dat) is included with the installed software.

A progress bar indicates the status of the firmware copying on the device, after which the actual update process begins. A counter appears on the display, which shows the number of data blocks remaining before the update process is completed. After the update is completed, a message appears on the display about the result of the action taken: in case of successful completion, the device will automatically restart with a new program, if an error occurs, it is necessary to reprogram the device.

Permitted Functions

Firmware for each of the four SD system models Plus are accessed using a special hardware key (which can be seen on the first screen image of the device display, shown in the *USER INTERFACE section* on page 14), which identifies the characteristics of the amplifier and its allowed functions. More specifically, these keys consist of a six-digit alphanumeric code:

- EP (Everted Potential) channels (46 for the Research model , 36 for the Flexi model and its two variants)
- The third digit indicates the presence (1) or absence (0) of a Bluetooth module (Bluetooth is not available on any of these models)
- The fourth character is a letter, and together with the first two digits of the code, identifies the device model (46 A = Research , 36 A = Flexi , 36 B = HighRate , 36 D = Clinic)
- The fifth character is a letter that indicates whether the device has a Micromed product label (A) or is distributed by a third party (C).
- The last character is the letter (E) which indicates the presence of two EP (Everted Potential) channels.

Thus, the four keys are associated with the four models as follows:

46OAAE	SD Plus Research
36OAAE	SD Plus Flexi
36OBAE	SD Plus Flexi HighRate
36ODAE	SD Plus Flexi Clinic

SD System Hardware Keys Plus

CHECKING THE DEVICE STATUS

SD System Amplifier Specifications PLUS can be controlled by specific tests using laboratory instruments and a data acquisition program. There is also a software package called **MICROMED DEVICE TEST (MICROMED DEVICE TEST)** . which includes a special program called **SD PLUS TEST (SD SYSTEM TEST PLUS)** , which automates the control process and generates a test report.

MANAGEMENT OF GENERAL FUNCTIONS

With respect to the function control described below, it is recommended to use connectors with 10 k Ω resistors, which are installed on the inputs of the amplifier. In addition, a voltage divider function generator must be used that generates the proper signals to feed into the amplifier.

It is necessary to have a function generator with an external trigger that allows regulation of the number of cycles (periods) of the signal.

RCA socket located on the interface between the personal computer and the amplifier, or in the window of the interface card. In the signal generated by the generator, it is necessary to set the number of cycles equal to five.

It is necessary to include a voltage divider with a ratio of 1:1000 in the system.

Measurements must be made using the SystemPlus software EVOLUTION .

MICROMED can supply an accessory kit that will allow you to make all the connections required for the tests described below .

CMRR (Common Mode Rejection Ratio) Test - EEG Channels


Figure 10 - Wiring diagram for CMRR tests

CMRR (Common Mode Rejection Ratio) test is performed by applying a symmetrical sinusoidal signal (make sure that the bias voltage in your device is zero), 50 Hz, 1.6 V amplitude for all EEG channels, including G 2 (negative input of both channels data acquisition) and G 1 (internal voltage reference). Set the EEG acquisition software to ALL mode CHANNELS (ALL CHANNELS), and set the sampling rate to 512 Hz.

Set filters to maximum bandwidth in high frequencies, 0 Hz in low frequencies, 1000 Hz capture window: 50 μ V/cm

Disable the notch filter and start the data collection process.

Based on the received data, it is possible to calculate the CMRR (Common Mode Rejection Ratio) in decibels using the formula:

$$CMRR = 20 \log_{10} \left(\frac{1,6}{\nu} \right)$$

where B is the signal amplitude from peak to peak in each channel.

This parameter must be equal to or greater than 125 dB.

The correspondence between some CMRR values and the measured amplitude is shown below.

μ V pk- pk	CMRR (when VIN=1.6V)	μ V pk-pk	CMRR (when VIN=1.6V)	μ V pk- pk	CMRR (when VIN=1.6V)
0.1	144	2	118	12	102
0.2	138	3	114	13	101
0.3	134.5	4	112	14	101
0.4	132	5	110	15	100
0.5	130	6	108	16	100
0.6	128.5	7	107	17	99
0.7	127	8	106	18	98
0.8	126	9	104	19	98
0.9	125	10	104	20	98
1	124	eleven	103	21	97

Table 6 - Correspondence between CMRR and measured voltage

CMRR Test (Common Mode Rejection Ratio) - Polygraphic Channels

CMRR (Common Mode Rejection Ratio) test is performed by applying a symmetrical sinusoidal signal (make sure that the bias voltage in your device is zero), a frequency of 50 Hz, an amplitude of 1.6 V for all 11 printing channels. Set the EEG acquisition program to C 21 P mode , and set the sampling rate of all 11 polygraphic channels to 1024 Hz.

Set filters to maximum bandwidth.

Output the amplitude of the signal using the video.

With this information, it is possible to calculate the CMRR in dB using the formula and values in Table 6:

This parameter must be greater than or equal to 125 dB.

Background noise test - EEG channels

Connect all inputs of EEG channel 21 and G 2 to G 1.

Set the sampling rate to 256 Hz, and preset the filters to a bandwidth of 0.3 - 30 Hz. Collect data and verify that the background noise is less than $0.9 \mu\text{V}$ pk-pk ($0.15 \mu\text{V}$ rms) in a 0.5 second time frame.

Background noise test - printing channels

Short all inputs (positive and negative) of 11 printing channels to G 1.

Set the sampling rate to 256 Hz, and preset the filters to a bandwidth of 0.3 - 30 Hz. Make sure the background noise is less than $0.9 \mu\text{V}$ pk-pk ($0.15 \mu\text{V}$ rms).

Gain Validity Test - EEG Channels

SDPLUS system

Divisor 1:1000

Signal generator


Figure 11 - Wiring diagram for gain tests

Set the signal generator to deliver a 5 Vp-p, 10 Hz square wave. Use a symmetrical voltage divider to feed the generated signal attenuated by a factor of 1000 ($5000 \mu\text{V}$ pk-pk 10 Hz is required) to the amplifier's polygraphic inputs. Connect the positive output of the divider to the positive inputs of the amplifier (channels 1 to 24), the negative output of the divider to the G 2 inputs, and to ground at G 1.

BQ software, set up a program that shows all EEG channels relative to G 2. In this program, the filters should be set to the largest bandwidth: 1024 Hz. Set the data collection function to "all channel" ("all channels") and a data acquisition window at $6400 \mu\text{V}$. Verify that the observed amplitude is $5000 \pm 500 \mu\text{V}$.

Gain Correctness Test - Polygraphic Channels

Set the signal generator to deliver a 5 Vp-p, 10 Hz square wave. Use a symmetrical voltage divider to feed the generated signal attenuated by a factor of 1000 ($5000 \mu\text{V}$ p-p 10 Hz should be obtained) to

the amplifier's polygraphic inputs (channels 25 - 40). connect the negative output of the divider to the negative inputs of the amplifier, and to ground at slot G 1. In the acquisition software, set the program to show all printing channels. In this program, the filters must be set to the largest bandwidth. Set the sampling rate to 128 Hz. Set the data collection function to "all channel" ("all channels") and a data acquisition window at 6400 μ V. Verify that the observed amplitude is $5000 \pm 500 \mu$ V.

Testing of EP channels (evoked potentials) with a large bandwidth

Set the function generator to supply a sinusoidal signal of 5 Vp-p, 10 Hz. Use a symmetrical voltage divider to apply a 1000 times attenuated signal to the EP 1 and EP 2 inputs. Verify that the displayed amplitude is $5000 \pm 500 \mu$ V.

Impedance check

Install 10 k Ω resistors between the 40 EEG/EP channel inputs and ground (G 1). Make sure the impedance checker reads the values correctly.

Check inputs EP 1 and EP 2 in the same way.

Automatic test software

Upon request, automatic test software for the **SD system amplifier family is available. PLUS .**

The software allows the user to perform background noise and gain correctness tests, input impedance measurements, and control functions of the built-in filter change switches and acquisition windows. In addition, appropriate connectors for the amplifier inputs are available to speed up testing. The program requires the use of a function generator.

ELECTROENCEPHALOGRAPHIC SYSTEM FOR NEUROPHYSIOLOGICAL STUDIES MICROMED , VERSION: MORPHEUS for PSG:

- Inspect the product, accessories, and recording cables for possible damage caused by transportation. Such problems should be checked and reported to your Micromed representative before using the device.
- Before cleaning the system, make sure it is turned off. Do not use abrasives or solvents of any kind.
- Clean the device with a "damp" cloth (not soaked) on the outside of the lid. Wipe the connection cables with a damp cloth. Replace electrodes after each use when necessary for hygienic reasons.
- If the device has been dropped from a height of less than 1 m, verify regular device functionality with a RESISTANCE CHECK test and calibration resistors before use on a patient. If the device is dropped to a greater height, Micromed personnel must check the device.
- Do not open the cover of the system or expose it to abnormal conditions (high temperature or excessive humidity).
- Avoid spilling any materials, in particular liquids, on the system components, and if this happens, do not use the system for any reason and contact your representative for the necessary assistance.
- There are no appliances inside that can be serviced by the user. Serviced by Micromed and authorized parties only. The warranty is void if the device is opened.

- Do not repair or replace any part of the devices arbitrarily. For any questions, please contact qualified personnel. In case of potentially dangerous events, check the device, even if it seems to be working well.
- Contact your local representative to arrange for a periodic maintenance review (at least every two years).
- Do not dispose of the device and used batteries with normal waste. Check your local electronic equipment disposal regulations.

1.2.21 The EEG system must ensure operability and electrical safety in case of contact with biological fluids on its surface.

The design must be made in such a way that it is impossible to touch live parts.

1.2.22 The design of the EEG system must be suitable for disinfection, ensure the flow of washing liquids; the number of grooves and protrusions that contribute to the accumulation of contaminants or biological fluids (sweat) should be minimal.

Disinfection of the surfaces of devices included in the system should be carried out by a chemical method, with a frequency of 1 time per day and after each use. As disinfectants in the territory of the Russian Federation, the manufacturer recommends using products registered in the register of the FGUZ "IMC Rosspotrebnadzor":

- 1) Dezol-Ch, with active ingredient - ethyl alcohol - 93.0 $\pm$ 3.0%;
- 2) "North", with the active ingredient - Ethyl alcohol - 93% at the rate of (93 $\pm$ 3)%, alkyldimethylbenzylammonium chloride - 0.07% at the rate of (0.07 $\pm$ 0.01)%;
- 3) "Anti-virus", with the active substance - Ethyl alcohol - 56%, alkyldimethylbenzylammonium chloride - 0.3% (HOUR);
- 4) "HELIOS", with the active substance - polyhexamethylene anidine hydrochloride - 0.1%, ethyl alcohol - 93%;
- 5) "DEZOL", with the active substance - olihexamethylene guanidine chloride (PHMG) - 0.08%, ethyl alcohol - 90%.

In the design of the EEG system, there should be no non-washable places, partitions, sharp narrowing of the cross section, and other inaccessible places.

1.2.23 The length of the lead cable from the connector for connecting to the EEG system to the electrodes must be at least 2.5 m.

1.2.24 The EEG system must comply with the requirements of this technical document if there is a constant voltage input on the corresponding tested channel $\pm(300 \text{ mV} \pm 10\%)$ between any lead electrodes.

1.2.25 The connector used to connect the patient's cable and the electrode lead connector should not have, when disconnected from the EEG system, any conductive parts connected to the patient's body that could come into contact with a flat conductive surface with a diameter of at least 100 mm.

Note - This requirement does not apply to the electrodes themselves .

1.2.26 Leakage current to the patient in the working ground circuit should not exceed 5 mA.

1.2.27. The software is pre-installed, the same for all types of system execution and must comply with the following characteristics:

Software safety class: Class B according to IEC 62304:2006.

Name of the software	Identification name of the software	Version number (identification number) of the software	Digital software identifier (checksum of the executable code)	Algorithm for Calculating Digital Software ID	Software date
System PLUS EVOLUTION - BRAIN QUICK	-	1.05.1	UDI : (01)0833928120 880(8012)1.05.1	-	06/15/2020

ATTENTION! The software options listed in the standard equipment of the systems (Appendix A) are pre-installed in all types of system versions without exception. However, they will remain inactive until the purchase of an access key for a particular software option.

Programming interface, may include:

Software option for registration and viewing:

- 2.1 Software option for recording and viewing EEG, 2-channel CFM version (on request)
- 2.2 Software option for recording and viewing EEG, 8-channel CFM version (on request)
- 2.3 Software option for recording and viewing EEG, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
- 2.4. Software option for recording and viewing digital video-EEG up to 128 channels, version: CLINIC (on request)
- 2.5. Software option for recording and viewing digital video-EEG up to 256 channels, CLINIC version (on request)
- 2.6 Software option for recording and viewing digital video-EEG up to 64 channels, version CLINIC (on request)
- 2.7 Software option for recording and viewing digital video-EEG, version CLINIC (on request)
- 2.8. Software option for recording and viewing digital video-EEG in Full HD format (on request)
- 2.9 Software option for recording and viewing digital video-EEG over the network, version CLINIC (on request)
- 2.10 Software option for recording and viewing digital video-EEG over a network of up to 128 channels, version CLINIC (on request)
- 2.11 Software option for recording and viewing digital video-EEG over a network up to 256 channels, CLINIC version (on request)
- 2.12 Software option for recording and viewing digital video-EEG over a network of up to 64 channels, version CLINIC (on request)
- 2.13 Software option for real-time logging and viewing over the network (as required)
- 2.14 Software option for EEG analysis in "on-line" and "off-line" modes, version: full / CFM / DSA. (on demand)
- 2.15 Software option for transmitting EEG data via TCP/IP (optional)
- 2.16 Software option "ERP -manager" (optional)
- 2.17 Software option "Management of outpatient registrar" (as required)
3. Software options for viewing:
 - 3.1. EEG software viewing option, version: CLINIC, STUDIO, CLINIC MULTICH (on request)
 - 3.2 Software option for viewing EEG, version: EASY (on request)
 - 3.3 Software option for viewing EEG in "off-line" mode, version: full, CFM, DSA (on request)
 - 3.4 Software option for connecting to databases of external applications (optional)
 - 3.5 EEG software option "Event identification" (as required)
 - 3.6 EEG software option "Analysis of high-frequency oscillations" (on request)
 - 3.7. Software option EEG "Brain Mapping", version: CLINIC, STUDIO (on request)
 - 3.8. EEG software option "Data exchange with MRI" (on request)
 - 3.9. Software EEG option "Interface HOLBERG SCORE" (optional)

- 3.11. Software EEG option "Averaging in "off-line" mode" (on demand)
- 3.12. Software option for viewing digital video-EEG, version: CLINIC, STUDIO (on request)
- 3.13. Software option for viewing digital video-EEG in Full HD format (on request)
- 3.14. Software option for viewing digital video-EEG in multichannel mode, version: CLINIC (on request)
- 3.15. Software option for advanced management of the archive server (optional)
- 3.16. Program statistics option (optional)
- 3.17. Software option for sleep program interface (optional)
- 3.18. Software option for background data transfer (optional)
- 3.19. Software option for exporting EEG data to the Viewer (optional)
- 3.20. Software option to enable HIS interface (optional)

4. Software options for sleep analysis:

- 4.1 Software option for analyzing imported curves in EDF format (optional)
- 4.2 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring SD LTM 32 EXPRESS (32 channels) (on request)
- 4.3 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring SD LTM 64 EXPRESS (64 channels) (optional)
- 4.4 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring 32 channels (on demand)
- 4.5 Software option for analyzing curves recorded on EEG Head Units for long-term monitoring 64 channels (on request)
- 4.6 Software option for the analysis of curves recorded on the EEG Head Units of the SD PLUS series (on request)

5. Program control options.

- 5.1 Software option "Centralized control of video cameras" (as required)
- 5.3 Software option "Event management" (as required)

ATTENTION! Information about the operation of each option and the operation of the software as a whole is reflected in the user manual.

Software signal reception and digitization

- Electroencephalography (EEG), long-term monitoring (LTM), EEG (CFM)
 - SD LTM 64 EXPRESS, SD LTM 32 EXPRESS
 - MORPHEUS
 - SD PLUS RESEARCH, SD PLUS FLEXI, SD PLUS FLEXI Clinic, SD PLUS FLEXI Highrate
- Evoked Potentials (EP)
 - SD PLUS RESEARCH, SD PLUS FLEXI, SD PLUS FLEXI Clinic, SD PLUS FLEXI Highrate
- Examination of the patient during sleep - polysomnography (PSG) and EEG
 - MORPHEUS
- Stimulants
 - Stimulator SD LTM STIM

Software data interfaces

BQ PCI EXPR LTM HALF/FULL: PCI bus interface for controlling amplifiers via cable; equipped with dual connection for Micromed devices and trigger input/output.

BQ USB: external interface connected to the USB port for controlling amplifiers via cable, with trigger in/out

BQ USB EXPRESS: external interface connected to the USB port for controlling amplifiers via cable, with trigger in/out and isolated serial port

BQ NET POE 2: network interface (USB converter/network based on BQ USB EXPRESS)

BQ NEL POE M: network interface (USB converter/BQ USB based network)

BOX TERMINAL EXPRESS: interface repeater, includes power supply adapter for medical equipment

User interface (monitor, keyboard, mouse)

LCD monitor 17" - 19" - 21", minimum resolution 1024x768

Keyboard

Mouse

PC/laptop/monoblock

Digital signal processing and storage software

Stationary PC or laptop or monoblock compliant with IEC 60950 and with minimum system requirements determined by the OS and installed version of Micromed software

Writer CD / DVD-drive (built-in PC / laptop / monoblock)

Video systems for full software operation

Standard digital video camera of the following models:

CAMERA IP HD PRO (article (serial number) Axis m5525)*

CAMERA IP HD OUTDOOR (article (serial number) Q 6075) *

CAMERA IP HD PLUS (article (serial number) V 5915) *

CAMERA IP COMPACT (article (serial number) M 1065- L)*

*Note: the articles are provided by the manufacturer's QMS, according to the manufacturer's QMS, the article (serial number) is indicated in the operational documentation, the article (serial number) is indicated on the nameplate of the device . However, the model is indicated on the marking of the camcorder on the body of the device.

Professional digital video camera with remote control and infrared illumination device

Standard graphics card or high resolution graphics card

Video switcher

Possibility and methods of integration with other medical devices .

The system can be used in conjunction with **Magstim magnetic stimulation devices** manufactured by The Magstim Company Limited (The Magstim Company Limited) , UK (not part of the system. Registration certificate No. RZN 2016/4594 dated 08/12/2016). To connect the system to stimulators, a CB MAG RAPID cable is provided for connecting the Micromed system to Magstim magnetic stimulation devices (D-Sub 26 pin (optional). The cable is included in the system.

1.3. Requirements for materials and components

1.3.1 The list of components, devices, instruments and accessories that make up the EEG system is presented in Appendix A.

1.3.2 The product must be made of the following materials: external parts of the components: GWPLAST 75 technopolymer, ABS brand TERLUX 2802 TR, buttons of the components (if any) silicone rubber brand FORSIL, cables, wires, extension cords, adapters: polyurethane brand BASF 60% , soft gauges, fixing tape: nylon brand Twill, PVC grade C-7058-M , brackets, chassis, holders, fasteners, trolleys. : AISI304 steel, PVC tongs C-7058-M , outer part of ABS sensors TERLUX 2802 TR.

Electrode caps: profiled silicone rubber brand FORSIL.

The electrodes are made of gilded plates or synthesized Ag / AgCl plates in a plastic sheath (BM polycarbonate).

Disposable ECG electrodes and have an adhesive (adhesive) layer: the non-woven component is a mixture consisting of 70% viscose and 30% HDPE (high pressure) impregnated with a hypoallergenic adhesive (adhesive composition (adhesive layer): a mixture of sodium carboxymethyl cellulose - from 50 to 60%, pectin - from 10 to 15% , gelatin - from 25 to 40%) in the amount of 30 g / sq.m. Strength of adhesion to the skin: It is provided with a soft adhesive layer of the adhesive surface, the adhesive surface gently adheres to the skin ensuring tight adhesion, while it should be noted that when the fastener is removed, the sticky layer does not cause pain or discomfort to the patient. the adhesive layer peels off gently and painlessly.

The force applied to remove the electrode is 2.5 N.

The connecting parts are only: electrodes, caps and sensors, tapes.

Descriptions of the material composition of each component system are given below.

The product must be made of the following materials:

outer parts of the component parts:

USB dongle for registration and viewing programs Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

USB dongle for viewers for 1/5/10/25 users Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

USB dongle for sleep analysis programs Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

CB TERMINAL 15 cable for connecting the interface board to the BOX TERMINAL EXPRESS interface module

Wall panel EXPRESS WALL CONNECTOR /12 with additional socket and IN-trigger and OUT-trigger for use in combination with BOX TERMINAL EXPRESS interface module (optional) and 12V voltage amplifier (optional) AISI304 steel grade

LCD display: 19"/22"/24"/27" (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS , FORSIL brand silicone rubber

EEG headunit for long-term monitoring SD LTM 32 EXPRESS (32 channels) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

Electrode box JB LTM 32 EXPRESS to EEG head unit for long-term monitoring SD LTM 32 EXPRESS (32 channels) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

BQ interface module USB PLUS to EEG headunit for long-term SD monitoring LTM PLUS (64 channels) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

EEG head unit for long-term SD monitoring LTM 32 PLUS (32 channels) GWPLAST 75 technopolymer, ABS grade TERLUX 2802 TR

16JB LTM 32 PLUS Electrode box to the EEG head unit for long-term SD monitoring LTM 32 PLUS (32 channels) GWPLAST 75 technopolymer, ABS grade TERLUX 2802 TR

Electrode box JBD SUB PLUS with Dual 37-Pin D - Sub Connector to EEG Head Units for SD Series Long-Term Monitoring LTM PLUS Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

JB TOPHAT PLUS electrode box for all channels with dual 37-pin D-Sub connector to EEG headunits for long-term monitoring SD series LTM PLUS Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

BOX TERMINAL MULTI interface module for connecting the EEG Head Unit for long-term SD monitoring LTM PLUS (up to 256 channels) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

BQ NET POE PLUS interface module for connecting EEG headunit for SD long-term monitoring LTM PLUS , includes SUPP wall mount BQ NET POE 2

CB LTM DC cable for connecting DC channels to EEG Head Units for long-term SD monitoring LTM GWPLAST 75 technopolymer, TERLUX 2802 TR ABS , BASF 60% polyurethane

EEG head unit for long-term monitoring SD LTM 64 EXPRESS (64 channels) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

Electrode box JB LTM 64 EXPRESS to EEG head unit for long-term monitoring SD LTM 64 EXPRESS (64 channels) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

SD LTM STIM Cortical Stimulator for EEG Head Units for Long-Term Monitoring SD Series LTM (64 channels) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS , FORSIL brand silicone rubber

Electrode box JB LTM 64 PLUS to EEG head unit for long-term monitoring SD LTM 64 PLUS (64 channels) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

EEG head unit for long-term monitoring SD LTM 64 PLUS (64 channels) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

EEG headunit SD PLUS CLINIC (29 channels) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

EEG headunit SD PLUS FLEXI (38 channels) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

EEG headunit SD PLUS HIGH RATE (29 channels) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

EEG SD PLUS RESEARCH headunit (48 channels) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

MORPHEUS EEG head unit for PSG Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

Electrode box JB MORPHEUS (24 channels) to EEG head unit MORPHEUS for PSG Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

BQ NET POE M interface module for connecting MORPHEUS EEG headunit for PSG, includes SUPP BQ NET POE 2 wall mount and BQ USB interface module GWPLAST 75 technopolymer, TERLUX brand ABS 2802 TR

BQ USB interface module to MORPHEUS EEG headunit for PSG (on request) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

BOX TERMINAL interface module for MORPHEUS EEG headunit for PSG, including FO PLUS circuit and SD power supply (on request) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

CB DC cable for connecting DC channels to MORPHEUS EEG headunit for PSG (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS , BASF 60% polyurethane

CB JB MORPHEUS cable for JB MORPHEUS electrode box polyurethane brand BASF 60% ARM GM MORPH bracket for MORPHEUS EEG head unit for PSG assembly (on request) AISI304 steel

BAG MORPH bag for MORPHEUS EEG head unit for PSG

Belt to the bag of the EEG head unit MORPHEUS for PSG

SUPP MORPH holder for MORPHEUS EEG head unit for PSG (on request) AISI304 steel

MORPH BUTTON patient button for use with MORPHEUS EEG headunit for PSG Silicone rubber brand FORSIL, ABS brand TERLUX 2802 TR

Oximetry sensor connection unit for use with MORPHEUS EEG headunit for PSG (on request) GWPLAST 75 technopolymer, TERLUX 2802 TR grade ABS

Electrode box JB LTM DSUB with double 37-pin D-Sub connector for EEG head units for long-term monitoring of the SD LTM EXPRESS series (on request) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

Electrode box JB LTM TOPHAT (universal) with sensor inputs for all channels (as required) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

LED photostimulator with bracket GWPLAST 75 technopolymer, TERLUX 2802 TR ABS , FORSIL silicone rubber

BQ NET POE 2 interface module for connecting SD LTM EXPRESS long-term EEG headunit, includes SUPP BQ NET POE 2 wall mount and BQ USB EXPRESS interface module (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

Wall panel EXPRESS WALL CONNECTOR / DE with additional socket and IN-trigger and OUT-trigger for use in combination with interface module BOX TERMINAL EXPRESS (on request) AISI304 steel

BQ NET POWER POE power supply for BQ NET POE series interface modules (on request) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

BQ USB EXPRESS interface module to EEG head unit for long-term monitoring SD LTM EXPRESS (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

BQ USB EXPRESS interface module to EEG head unit for long-term monitoring SD LTM EXPRESS (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

BQ USB PLUS interface module to EEG headunit for long-term monitoring SD LTM EXPRESS (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

USB adapter to serial port ADAPT USB SER (as required) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

Isolated USB to Serial Adapter ADAPT USB SER ISO (as required) Technopolymer GWPLAST 75, ABS grade TERLUX 2802 TR

CB LTM DC cable for connecting DC channels to EEG Head Units for long-term SD monitoring LTM polyurethane brand BASF 60%

CB MAG RAPID cable for connecting the Micromed system to Magstim magnetic stimulation devices (D-Sub 26 pin (optional) BASF brand polyurethane 60%

Trigger cable CB RCRC4N for stimulator (optional) BASF brand polyurethane 60%

Trigger cable CB RCRC4R for stimulator (optional) BASF brand polyurethane 60%

Trigger cable CB RCRC1R for stimulator (optional) BASF brand polyurethane 60%

Pair of CB RCRC5 PAIR Trigger Cables for Stimulator (as required) BASF brand polyurethane 60%

Cable CB SER RCA EXT for synchronization with external stimulators (on demand) polyurethane brand BASF 60%

Fiber optic cable 10m/15m/20m/30m/40m/50m/100m (as needed) BASF brand polyurethane 60%

HBC EXPRESS cable for connecting an interface board or interface module to EEG Head Units (on request) BASF grade polyurethane 60%

HBC INT LTM R cable for connecting an interface board or interface module to SD LTM EEG Head Units for long-term monitoring (as required) BASF grade polyurethane 60%

Power supply LED POWER for photostimulator LED Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

TERM LTM Quick Disconnect Cable for EEG Head Units polyurethane brand BASF 60%

PRO (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS , FORSIL brand silicone rubber

CAMERA HD Full HD (1920 x 1080) with PTZ control and automatic IR switching (on demand) (on demand) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS, FORSIL brand silicone rubber

CAMERA IP HD OUTDOOR (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS , FORSIL silicone rubber

CAMERA IP HD PLUS video camera (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS , FORSIL brand silicone rubber

including microphone (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS , FORSIL brand silicone rubber

CAMERA IP COMPACT (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS , FORSIL silicone rubber

HDMI EXTENDER device for transmitting HDMI signal over Cat 5 e /6 cable, including receiver and transmitter (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

IP MIC microphone for IP video cameras (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS , FORSIL brand silicone rubber

IP POWER power adapter for IP video cameras (on request) BASF 60% polyurethane,

Composite switcher VMIXER PRO with 2 inputs (optional) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

Professional switcher VMIXER PRO with 5 HDMI inputs (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

Bracket ARM GM LTM for EEG head units for long-term monitoring of the SD LTM series, assembled (on request)

Vertical bracket ARM LCD DOUBLE for two LCD displays for fixing on a trolley (as required) AISI304 steel

Horizontal bracket ARM LCD HORIZ for two LCD displays for fixing on the trolley (on demand) steel grade AISI304

Vertical bracket ARM LCD SW for LCD displays for fixing on the trolley (on demand) steel grade AISI304

Adjustable BAG BELT HP to secure the headunit bag to the patient

EXPRESS WALL CONNECTOR NT /12 wall panel with additional socket for use in combination with BOX TERMINAL interface module (optional) EXPRESS and 12V voltage amplifier (optional) AISI304 steel

BAG LTM HP bag for EEG head unit for long-term monitoring SD LTM series (on request)

Rack with base BASE STAND (on request) AISI304 steel

Clamp for fastening the head unit BQ CLAMP (on demand) steel grade AISI304

Swivel mount CAM JOIN for a video camera on a stand (optional)

CAM STAND for mounting a video camera or head unit on a trolley (optional) AISI304 steel

Fixing FIX GM for the bracket on the trolley (on demand) steel grade AISI304

Fixing FIX GM TW for bracket on the table (optional) AISI304 steel

JOIN GM connection for mounting the head unit holder on the bracket (optional) AISI304 steel

Screw LTM SCREW for attaching the electrode box to the head unit (optional) AISI304 steel

LED STAND PHOTIC wheelbase for LED photostimulator (on request) ABS brand TERLUX 2802 TR

Bracket LED SUPP PHOTIC for fixing the LED photostimulator (on demand) AISI304 steel

Holder STAND CAM SW video camera assembly (on demand) AISI304 steel

Rack STAND POLE NEW on wheels with weighting agent (on demand) steel grade AISI304

Wall mount SUPP BQ NET POE 2 wall mount for BQ NET POE 2 interface modules BQ NET POE 2 / BQ NET POE M (optional) (optional) (optional) AISI304 grade steel

Holder SUPP LTM for EEG Head Units for long-term monitoring of the SD LTM series (on request) steel grade AISI304

VCLAMP mount for video camera or IR illumination (optional) AISI304 steel

Tripod stand VSTAND for video camera or IR illumination (on request) AISI304 steel

VSUPP wall bracket for video camera (optional) AISI304 steel

Metal support VSUPP IP for attaching an IP video camera to a swivel mount (as required) steel grade AISI304

Metal support VSUPP IP IR for attaching the IP video camera to the rack (optional) AISI304 steel

Wall mount VSUPP IP HD PLUS BLACK for video camera, black color (on request) AISI304 steel

Wall mount VSUPP IP HD PLUS WHITE for video camera, white color (on request) AISI304 steel

Adapter ADAPT BNC RCA (on request) polyurethane brand BASF 60%

Adapter ADAPT DVI HDMI (on request) polyurethane brand BASF 60%

ANY WHERE 2 device for USB over LAN connection, with 2 USB ports (optional) (optional) ABS brand TERLUX 2802 TR

ANY WHERE 5 device for USB connection via LAN, with 5 USB ports (as required) ABS brand TERLUX 2802 TR

USB joystick CAM JOY STICK for camera control (optional) FORSIL silicone rubber, GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

24-key CAM PAD numeric keypad with USB connection to control cameras (optional) ABS brand TERLUX 2802 TR

Numeric keypad CAM PAD 48 by 48 keys with USB connection to control cameras (on demand) ABS brand TERLUX 2802 TR

Connecting cable CB BBNB1 (BNC - BNC) (optional) BASF brand polyurethane 60%

Connecting cable CB BNRC1 (BNC - RCA (Red) (as required) BASF brand polyurethane 60%

Connecting cable CB BNRC5 (BNC - RCA (Black) (on demand) polyurethane brand BASF 60%

CB DVI monitor cable (optional) BASF brand polyurethane 60%

Cable HDMI - HDMI MINI for video signal transmission, length 2m/5m/10m (as needed) polyurethane brand BASF 60%

HDMI to HDMI cable for high definition video transmission, length 2m/5m/15m (optional) BASF brand polyurethane 60%

Audio cable CB JACK JMINI with connectors 6.3 mm - MiniJack 3.5 mm (optional) BASF brand polyurethane 60%

Audio cable CB JACK RCRC with 6.3 mm plugs - double RCA (optional) BASF brand polyurethane 60%

Microphone cable CB MICROPHONE (optional) BASF brand polyurethane 60%

CB MON monitor cable (D-Sub 15 pin - D-Sub 15 pin) (optional) BASF brand polyurethane 60%

Two-core audio cable CB RCRC JMINI, RCA male-MiniJack 3.5 mm (as required) BASF brand polyurethane 60%

Two-wire audio cable CB RCRC PAIR, RCA male-double RCA (on request) polyurethane brand BASF 60%

CB VIDEO video cable as required) BASF brand polyurethane 60%

CB IP HD PLUS 15 cable for controlling IP HD PLUS cameras (optional) BASF grade polyurethane 60%

CB IP HD PLUS DC 5 cable for connecting IP HD PLUS cameras to each other (on request) BASF grade polyurethane 60%

CB IP HD PLUS EXTENDER adapter for IP HD PLUS signal transmission over Cat5e/6 cable, pair (optional) BASF brand polyurethane 60%

TERLUX 2802 TR ABS self-powered HDMI BOOSTER amplifier (on request)

Trolley CARSTAR mobile (on demand) ABS brand TERLUX 2802 TR

Bracket FIX TOR for attaching the isolation transformer to the trolley (optional) AISI304 steel

Mouse pad small/regular/large (on demand) PVC grade C-7058-M

Trolley CAR PLUS mobile, adjustable in height (according to needs) ABS brand TERLUX 2802 TR

FIX CART plates for fixing trolley supports (on request) AISI304 steel

Power cable set KIT CB POWER for trolley (on request) BASF grade polyurethane 60%

Pull-out keyboard shelf SUPP KEYB PLUS for trolley (on request) ABS brand TERLUX 2802 TR

Computer mount SUPP PC PLUS on a trolley (optional) AISI304 steel

Small computer mount SUPP PC PLUS SFF on trolley (optional) AISI304 steel

Horizontal bar BAR POLE NEW for mounting the bracket on the trolley (on demand) steel grade AISI304

Trolley CAR TORO with shelves, standard size (on request) AISI304 steel

Trolley CAR TORO SMALL with shelves, small size (on request) AISI304 steel

Rear cover COVER TORO for standard CARTORO trolley (optional) AISI304 steel

Back cover COVER TORO SMALL for small trolley CAR TORO SMALL (on request) steel grade AISI304

Isolation transformer TOR 500, Pmax 500 VA, black (on request) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

Isolation transformer TOR 500 W, Pmax 500 VA, white (on request) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

Isolation transformer TR 50, 100 VA (on request) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

Network isolator NET ISOLAT for networks 10/100 Mbit, with RJ45 connectors (on request) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

Network cable NET RJ45 with RJ45 connector, length 3m/5m/10m (on demand) polyurethane brand BASF 60%

NET SWITCH 4 4-port network switch (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

NET SWITCH 8 8-port network switch (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

Professional network switch NET SWITCH PRO (on demand) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

LTM BUTTON EXT extension for LTM BUTTON TP patient button, length 120 cm (as required) BASF brand polyurethane 60%

Patient Button LTM BUTTON TP for use with EEG Head Units for long-term monitoring SD LTM silicone rubber brand FORSIL, ABS brand TERLUX 2802 TR

Pliers TP tongs for removing electrodes from the electrode box (on demand) PVC grade C-7058-M

USB extension cable USB EXT CB (optional) BASF 60% polyurethane

USB EXTENDER device for transmitting a USB signal over a Cat5e/6 cable, including a receiver and transmitter (optional) ABS brand TERLUX 2802 TR

BQ CONTACT module for contact closure control (on request) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

Patient button with 5 m cable and RCA connector (optional) TERLUX 2802 TR ABS , FORSIL silicone rubber, BASF 60% polyurethane

BUTTON BOX for patient button control (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

Cables for BUTTON BOX (on request) BASF brand polyurethane 60%

Illuminated patient button with 2 m cable and LEMO connector (optional) TERLUX brand ABS 2802 TR, FORSIL brand silicone rubber, BASF brand polyurethane 60%

Power supply for BUTTON BOX (optional) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

Cable for connection unit for oximetric sensors. Length 3 m (as required) BASF brand polyurethane 60%

Oximetry Sensor Connection Unit for use with EEG Head Units for long-term monitoring of the SD LTM series (on request) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

Oximetry sensor "EarClip", adult (on demand) nylon brand Twill, PVC brand C-7058-M, outer part of ABS sensors brand TERLUX 2802 TR.

FingerClip oximetry sensor, adult (optional) nylon brand Twill, PVC brand C-7058-M, outer part of ABS sensors brand TERLUX 2802 TR.

Oximetric sensor "FingerClip", children's (on demand) nylon brand Twill, PVC brand C-7058-M, outer part of ABS sensors brand TERLUX 2802 TR.

Fixing tape SD OXY RA, adult size (on request) nylon brand Twill, PVC grade C-7058-M

Fixation tape SD OXY RC, child size (as required) brand nylon Twill, PVC grade C-7058-M

Flexible oximetry sensor, adult size (as required) nylon grade Twill, PVC brand C-7058-M,

brand TERLUX 2802 TR.

Flexible oximetry sensor, child size (as required) nylon grade Twill, PVC brand C-7058-M, outer part of ABS sensors brand TERLUX 2802 TR.

Large reusable soft oximeter sensor (digital readout 12.5 - 25.5 mm) (optional) nylon grade Twill, PVC brand C-7058-M, outer part of ABS sensors brand TERLUX 2802 TR.

Medium size reusable soft oximeter probe (digital readout 10 - 19 mm) (optional) Grade nylon Twill, PVC brand C-7058-M, outer part of ABS sensors brand TERLUX 2802 TR.

Reusable soft oximetry sensor small size (digital readout 7.5 - 12.5 mm) (optional) nylon grade Twill, PVC brand C-7058-M, outer part of ABS sensors brand TERLUX 2802 TR.

COMP DVI RX KVM extender compact, for DVI + USB HDI reception (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

KVM Extender COMP DVI RX AUD compact, for receiving DVI + USB HDI + analog audio signal (on request) GWPLAST 75 technopolymer, ABS brand TERLUX 2802 TR

KVM Extender COMP DVI TX compact, for DVI + USB HDI transmission (on demand) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

COMP DVI TX AUD compact KVM extender for DVI + USB HDI + analog audio transmission (on request) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

KVM DVI module RX for DVI + USB reception HID (on request) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

KVM DVI module RX AUD dual for DVI + USB reception HID + analog audio signal (as required) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

KVM DVI module TX for DVI + USB transmission HID (on demand)

KVM DVI module TX AUD dual for DVI + USB transmission HID + analog audio signal (as required) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

KVM DVI module TX AUD DCAT dual for DVI + USB transmission HID + analogue audio signal with two-channel transmission (on demand) GWPLAST 75 technopolymer, ABS brand TERLUX 2802 TR

KVM DVI module TX DCAT to DVI + USB transmission HID ,dual transmission (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

KIT KVM KVM module including transmitter and receiver for keyboard, mouse and monitor. (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

KVM module KIT KVM DUAL, including transmitter and receiver, for keyboard, mouse and two monitors (on request) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

KVM Chassis KVM BODY 2 for 2 DVI modules (on request) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

KVM Chassis KVM BODY 4 for 4 DVI modules (on request) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

KVM Chassis KVM BODY 4 R for 4 DVI modules , including rack mounting bracket (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

KVM Chassis KVM BODY 6 for 6 DVI modules (on request) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

143
KVM Chassis KVM BODY 6 R for 6 DVI modules including rack mounting bracket (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

IP KVM module, KVM IP RX, receiver (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

KVM IP module, KVM IP TX, transmitter (optional) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

KVM matrix KVM SWITCH for 8/16/32/48/64/80 ports (on demand) Technopolymer GWPLAST 75, ABS brand TERLUX 2802 TR

Additional mount SUPP TOPHAT for Electrode box JB LTM TOPHAT (on request) steel grade AISI304

EXPRESS WALL CONNECTOR NT /DE wall panel with additional socket for use in combination with the BOX TERMINAL EXPRESS interface module (on request) AISI304 steel

HBC BT/RJ45 extension adapter for BOX TERMINAL EXPRESS interface module (module side)(as required) BASF grade polyurethane 60%

HBC HT/RJ45 extension adapter for BOX TERMINAL EXPRESS interface module (PC side) (on request) BASF grade polyurethane 60%

External battery charger EXTERNAL CHARGER (on request) GWPLAST 75 technopolymer, TERLUX 2802 TR ABS

Cable HBC EXTERNAL BATTERY for connecting an external battery (on request) polyurethane brand BASF 60%

Extension cable LED EXT for LED photostimulator (on request) polyurethane brand BASF 60%

Quick disconnect cable TERM EXPRESS to the EEG Head Unit for long-term monitoring SD LTM EXPRESS in a configuration of up to 64 channels (on demand) polyurethane brand BASF 60%

Electrode cap (optional) profiled silicone rubber brand FORSIL:

- MRINEO19 X

- MRINEO31 X

- MRINEO61 X

Cable for electrodes EL CFM (as required) polyurethane brand BASF 60% :

- CB CFM050

- CB CFM150

Electrodes with cable for LTM (on request) BASF 60% polyurethane cables, electrodes made of gold-plated plates or synthesized Ag / AgCl plates in a plastic sheath (BM polycarbonate):

- CB EL TPCO 19

- CB EL TPCO 31

Electrode (plate) EL CFM for 2-channel EEG (on request) electrodes are made of gold-plated plates or synthesized Ag / AgCl plates in a plastic sheath (BM polycarbonate)

Disc electrode with cable (optional) electrodes are made of gold-plated plates or synthesized Ag / AgCl plates in a plastic sheath (BM polycarbonate):

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Pediatric disk electrode with cable (optional) electrodes are made of gold-plated plates or synthesized Ag / AgCl plates in a plastic sheath (BM polycarbonate):

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Ear electrode (clip) EL TPEAR with cable (optional) electrodes are made of gold-plated plates or synthesized Ag / AgCl plates in a plastic sheath (BM polycarbonate)

EL TPLINK cable for connecting two ear electrodes EL TPEAR (on demand) polyurethane brand BASF 60%

Connecting cable for LTM / LTM 256 (optional) polyurethane brand BASF 60%

CB PUD150 cable for EL PUD electrodes (on request) BASF brand polyurethane 60%

Disposable electrodes EL PUD (on demand) outer PVC grade C-7058-M,

the non-woven component is a mixture consisting of 70% viscose and 30% HDPE (high pressure) impregnated with a hypoallergenic adhesive (adhesive composition (adhesive layer): a mixture of sodium carboxymethyl cellulose - from 50 to 60%, pectin - from 10 to 15% , gelatin - from 25 to 40%) in the amount of 30 g / sq.m. Strength of adhesion to the skin: It is provided with a soft adhesive layer of the adhesive surface, the adhesive surface gently adheres to the skin ensuring tight adhesion, while it should be noted that when the fastener is removed, the sticky layer does not cause pain or discomfort to the patient, the adhesive layer peels off gently and painlessly.

The force applied to remove the electrode is 2.5 N.

Ag / AgCl plates in a plastic sheath (BM polycarbonate):

- EMG MBOT2530
- EMG MBOT2540
- EMG MBOT3745
- EMG MBOT5050
- EMG MBOT7570

Snoring sensor (optional) nylon grade Twill, PVC brand C-7058-M :

- EPM 1200
- EPM 1250

EPM 1312 respiration sensor kit, includes sensor (EPM 1387) and 4 tapes of various sizes (EPM 1340 + EPM 1341 + EPM 1342 + EPM 1343) (optional) Grade nylon Twill, PVC grade

145

Breath sensor pediatric EPM 1390 (on request) nylon brand Twill, PVC grade C-7058-M
Air flow sensor pediatric EPM 1400 (on demand) nylon brand Twill, PVC grade C-7058-M
GM / HM Disposable Breath Sensor Set , includes: EPM 1415 cable and 10 sensors (Adult 1426 EPMs / Medium 1427 EPMs / Infant 1428 EPMs) (optional) Nylon Grade Twill, PVC grade C-7058-M

Air flow sensor EPM 1450 (optional) nylon grade Twill, PVC grade C-7058-M

Set of air flow and pressure sensors EPM 1475 (optional) nylon grade Twill, PVC grade C-7058-M

Air pressure sensor and snoring EPM 14830 (on request) nylon grade Twill, PVC grade C-7058-M

Body position sensor EPM 1563 (optional) nylon grade Twill, PVC grade C-7058-M

Motion sensor kit EPM 1715 (including piezoelectric transducer (EPM 1770) and pairs of ankle straps (EPMS 1730) (optional) Grade nylon Twill, PVC grade C-7058-M

Chest force transducer kit including interface cable and 4 chest straps (optional) Nylon grade Twill, PVC brand C-7058-M , polyurethane brand BASF 60% :

- EPM 9102 A for adults
- EPM 9102 C pediatric

Abdominal force transducer kit including interface cable and 4 abdominal bands (optional) nylon grade Twill, PVC brand C-7058-M , polyurethane brand BASF 60% :

- EPM 9101 A for adults
- EPM 9101 C pediatric

Set of abdominal and thoracic force transducers, including abdominal and thoracic interface cables, and 4 straps (optional) nylon grade Twill, PVC brand C-7058-M , polyurethane brand BASF 60% :

- EPM 914x A for adults
- EPM 914x C pediatric
- EPM 915x A for adults
- EPM 915x C pediatric

Breath sensor tape (optional) nylon grade Twill :

- EPMs 1340
- EPMs 1341
- EPMs 1342
- EPMs 1343

Piezocrystalline breath sensor EPMs 1387 (as required) nylon grade Twill, PVC brand C-7058-M , ABS brand TERLUX 2802 TR

Cable for air flow sensor single channel EPMs 1415 (on request) nylon brand Twill, PVC brand C-7058-M , ABS brand TERLUX 2802 TR

Thermal breathing sensor (on demand) ABS brand TERLUX 2802 TR :

- EPMs 1426 adult size
- EPMs 1427 medium size
- EPMs 1428 pediatric

EPMs 1450-B interface for EPM 1450 air flow sensor (optional) TERLUX brand ABS 2802 TR

S transducer for EPM 1450 airflow sensor (on request) ABS brand TERLUX 2802 TR

Combined thermistor pressure and air flow sensor EPMs 1475 ABS brand TERLUX 2802 TR
CO 2 cannula (on demand) PVC grade C-7058-M :

- EPMs 14802, nasal and oral, male adult
- EPMs 14803, male for newborns
- EPMs 14804, male pediatric
- EPMs 14805, male adult

Piezoelectric pressure transducer for air flow and snoring EPMs 14830 (on request) ABS brand TERLUX 2802 TR

Elastic band EPMs 1550 for body position sensor EPM 1520 (optional) nylon brand Twill

EPMs 1563 transducer for body position sensor EPM 1563 (optional) ABS brand TERLUX 2802 TR

146

Cannula ThermoCan (on request) PVC grade C-7058-M :

- EPMs 15802

- EPMs 15807

Pair of elastic ankle bands EPMs 1730 (optional) nylon grade Twill

Piezoelectric transducer EPMs 1770 for motion sensor EPM 1715 (optional) ABS brand TERLUX 2802 TR

Pair of induction bands (as needed) nylon grade Twill :

- 9002-L120 EPMs

- 9002-L133 EPMs

- 9002-L40 EPMs

- 9002-L60 EPMs

- 9002-L90 EPMs

Adjustable induction band EPMs 9010 (as required) nylon grade Twill

Induction interface cable (as required) BASF brand polyurethane 60% :

- EPMs 9101 for the abdomen

- EPMs 9102 for the thoracic region

ECG electrode disposable ECG EL (on demand) external PVC brand S-7058-M,

the non-woven component is a mixture consisting of 70% viscose and 30% HDPE (high pressure) impregnated with a hypoallergenic adhesive (adhesive composition (adhesive layer): a mixture of sodium carboxymethyl cellulose - from 50 to 60%, pectin - from 10 to 15% , gelatin - from 25 to 40%) in the amount of 30 g / sq.m. Strength of adhesion to the skin: It is provided with a soft adhesive layer of the adhesive surface, the adhesive surface gently adheres to the skin ensuring tight adhesion, while it should be noted that when the fastener is removed, the sticky layer does not cause pain or discomfort to the patient, the adhesive layer peels off gently and painlessly.

The force applied to remove the electrode is 2.5 N.

ECG electrode cable (optional) BASF 60% polyurethane :

-ECG EL TPCLB black

-ECG EL TPCLR red

2. Requirements for environmental protection

2.1 EEG systems under normal operating conditions are not a source of environmental pollution.

During their manufacture, waste that is hazardous to the environment is not generated.

The materials used are recyclable.

2.2 The main type of possible hazardous impact on the environment is pollution of the atmospheric air of populated areas, soils and waters as a result of

- emergency leaks (placers) of production materials;

- unorganized incineration and disposal of waste on the territory of the manufacturing enterprise

the teacher or outside it;

- their arbitrary dumping in places not intended for this purpose.

2.3 EEG systems and materials used in their manufacture should not pose a danger to human life, health and the environment both during operation and after its completion.

5. Marking

5.1.1 EEG system must have a general marking that includes the following data:

- name of the manufacturer;
- address of the manufacturer
- details of the authorized/official representative (name, address, telephone number)
- full name of the EEG system;
- supply voltage;
- service symbols;
- serial number;
- degree and class of protection against electric shock;
- date of manufacture
- refer to the user manual
- **symbol "keep away from moisture"**
- **symbol "carefully fragile"**
- **mains voltage and frequency**
- **protection class and type of working part**
- **recycling information**
- life time
- number and date of registration certificate

5.1.2 Transport marking - with handling signs: "Fragile. Caution", "Top", and "Keep dry".

EEG head units, electrode boxes must be marked with the following data:

- name of the manufacturer;
- address of the manufacturer
- details of the authorized/official representative (name, address, telephone number)
- full name of the block;
- supply voltage;
- service symbols;
- serial number;
- degree and class of protection against electric shock;
- date of manufacture
- refer to the user manual
- **symbol "keep away from moisture"**
- **symbol "carefully fragile"**
- **mains voltage and frequency**
- **protection class and type of working part**
- **recycling information**
- **brief information about the battery**
- life time

The photostimulator must be labeled with the following data:

- name of the manufacturer;
- address of the manufacturer
- details of the authorized/official representative (name, address, telephone number)
- full name of the block;
- supply voltage;
- service symbols;
- serial number;
- degree and class of protection against electric shock;
- date of manufacture
- refer to the user manual
- **symbol "keep away from moisture"**
- **symbol "carefully fragile"**
- **mains voltage and frequency**
- **protection class and type of working part**
- **recycling information**
- **power**

- **lamp supply voltage**
- **flash period**
- **flash frequency**
- **life time**

The video camera must be marked with the following data:

- name of the manufacturer;
- address of the manufacturer
- details of the authorized/official representative (name, address, telephone number)
- full name of the block;
- supply voltage;
- service symbols;
- serial number;
- degree and class of protection against electric shock;
- date of manufacture
- refer to the user manual
- life time
- **symbol "keep away from moisture"**
- **symbol "carefully fragile"**
- **mains voltage and frequency**
- **power**
- **recycling information**

The video switcher must be marked with the following data:

- name of the manufacturer;
- address of the manufacturer
- details of the authorized/official representative (name, address, telephone number)
- full name of the block;
- supply voltage;
- service symbols;
- serial number;
- degree and class of protection against electric shock;
- date of manufacture
- refer to the user manual
- **symbol "keep away from moisture"**
- **symbol "carefully fragile"**
- **mains voltage and frequency**
- **recycling information**
- **life time**

The isolation transformer must be marked with the following data:

- name of the manufacturer;
- address of the manufacturer
- details of the authorized/official representative (name, address, telephone number)
- full name of the block;
- supply voltage;
- service symbols;
- serial number;
- degree and class of protection against electric shock;
- date of manufacture
- refer to the user manual
- **symbol "keep away from moisture"**
- **symbol "carefully fragile"**
- **mains voltage and frequency**
- **power**
- **recycling information**
- **maximum power consumption**
- **network connection**
- **fuse information**

- **connector data**
- **thermal protection data**
- **life time**

The network isolator must be marked with the following data:

- name of the manufacturer;
- address of the manufacturer
- details of the authorized/official representative (name, address, telephone number)
- full name of the block;
- supply voltage;
- service symbols;
- serial number;
- degree and class of protection against electric shock;
- date of manufacture
- refer to the user manual
- **symbol "keep away from moisture"**
- **symbol "carefully fragile"**
- **mains voltage and frequency**
- **recycling information**
- **dielectric strength data**
- **protection class and type of working part**
- **life time**

The network switch must be labeled with the following information:

- name of the manufacturer;
- address of the manufacturer
- details of the authorized/official representative (name, address, telephone number)
- full name of the block;
- supply voltage;
- service symbols;
- serial number;
- degree and class of protection against electric shock;
- date of manufacture
- refer to the user manual
- **symbol "keep away from moisture"**
- **symbol "carefully fragile"**
- **mains voltage and frequency**
- **recycling information**
- **power**
- **class and type of working part**
- **life time**

The contact closure control module must be marked with the following data:

- name of the manufacturer;
- address of the manufacturer
- details of the authorized/official representative (name, address, telephone number)
- full name of the block;
- supply voltage;
- service symbols;
- serial number;
- degree and class of protection against electric shock;
- date of manufacture
- refer to the user manual
- **symbol "keep away from moisture"**
- **symbol "carefully fragile"**
- **mains voltage and frequency**
- **recycling information**
- **current value**
- **surge protection**

- **pulse output frequency**
- **breakdown voltage**
- **life time**

The block for connecting oximetric sensors must be marked with the following data:

- name of the manufacturer;
- address of the manufacturer
- details of the authorized/official representative (name, address, telephone number)
- full name of the block;
- supply voltage;
- service symbols;
- serial number;
- degree and class of protection against electric shock;
- date of manufacture
- refer to the user manual
- **symbol "keep away from moisture"**
- **symbol "carefully fragile"**
- **mains voltage and frequency**
- **recycling information**
- **protection class and type of working part**
- **life time**

The KVM module must be marked with the following data:

- name of the manufacturer;
- address of the manufacturer
- details of the authorized/official representative (name, address, telephone number)
- full name of the block;
- supply voltage;
- service symbols;
- serial number;
- degree and class of protection against electric shock;
- date of manufacture
- refer to the user manual
- **symbol "keep away from moisture"**
- **symbol "carefully fragile"**
- **mains voltage and frequency**
- **recycling information**
- **life time**

The cortical stimulator must be labeled with the following information:

- name of the manufacturer;
- address of the manufacturer
- details of the authorized/official representative (name, address, telephone number)
- full name of the block;
- supply voltage;
- service symbols;
- serial number;
- degree and class of protection against electric shock;
- date of manufacture
- refer to the user manual
- **symbol "keep away from moisture"**
- **symbol "carefully fragile"**
- **mains voltage and frequency**
- **recycling information**
- **maximum output voltage**
- **isolation**
- **nutrition**
- **maximum stimulus intensity**
- **maximum operating voltage**

- 151
- maximum charge per phase
 - class/type of working part
 - life time

The label affixed to the packaging of electrodes, cables, sensors should contain the following information:

- name of the manufacturer;
- address of the manufacturer
- month, year of manufacture;
- serial number;
- best before date

6. Packing

6.1 The equipment of the EEG system is delivered complete packed in polyethylene film and in boxes placed in a shipping container - corrugated cardboard boxes.

The container must be closed and sealed with tape.

6.2 The gross weight of the shipping container should be no more than 33 kg.

6.3 Shipping documentation for EEG systems is supplied with them in a waterproof plastic bag.

The packing list must include:

- name of the manufacturer;
- EEG system name
- factory (serial) number;
- the number of products in the package;
- batch or order number;
- date of packing;
- completeness.

7. Repair details

In case of post-warranty and out-of-warranty repair of the product, please contact: Infomed-Neuro Limited Liability Company. jur. address: 125445, Moscow, Leningradskoe highway, 69, building 1, room 1, floor 1, room 50. phone 8 (495) 645-47-00. kitayeva@im-neuro.ru.

8. Safety requirements

8.1 EEG systems are safe when used in accordance with the rules of operational documentation, they are not a source of dangerous and harmful factors.

8.2 EEG systems must comply with the design safety and ergonomic standards of EN 614-1:1995, ISO 26800:2011, ISO/TR 16142:2006, EN 60601-2-26:2003, IEC 60601-1.

8.3 Electrical safety requirements are specified in the section Technical characteristics and Technical specifications.

9. Storage and transportation

9.1 Requirements for transportation and storage:

152

Transportation of EEG systems is carried out by all types of covered transport in accordance with the rules for the carriage of goods in force for this type of transport.

9.2 The conditions for the transportation of EEG systems in terms of the impact of climatic factors must comply with: temperature from -50 to +50 °C relative humidity 75-100%.

9.3 Storage conditions - temperature from -50 to +40 °C relative humidity 75-98%.

9.4 EEG systems are stored in conditions that exclude exposure to direct sunlight, precipitation and aggressive environments.

10. Instructions for use

10.1 EEG systems should be used for purposes strictly in accordance with the manufacturer's instructions.

Attention! When working with the system and its components, as well as when assembling the system, be sure to: wear disposable gloves to avoid personnel contact with the system and the transfer of infectious agents to external parts of the system. You also need to be wearing a mask.

The patient is also required to wear disposable gloves (it is allowed to remove gloves only when working with sensors that are attached to the finger). This requirement was introduced in connection with the control of the situation on COVID -19.

The contact parts are only: electrodes, caps and sensors.

10.2 The room where the EEG system is supposed to be used must be equipped with an electrical network with grounding.

10.3 When examining patients, the doctor must comply with the Code of Medical Ethics.

10.4 Operational documentation of EEG systems should contain information on indications and contraindications for use .

EEG systems are useful for studying a range of neurological diseases, evaluating psychiatric disorders, and helping locate tumors or lesions on or near the surface of the brain.

10.5 Main functionality:

- continuous monitoring and registration of subelectrode resistances and potentials;
- all types of mathematical processing - spectral, autocorrelation and cross-correlation analysis, coherence function with topographic mapping, mathematical statistics;
- analysis of the functional asymmetry of the brain with an assessment of the depth of interhemispheric and intrahemispheric interactions;
- automatic search for manifestations of paroxysmal activity during long-term studies for the diagnosis of epilepsy;
- automatic generation of the description and classification of the electroencephalogram with the possibility of editing;
- Maintaining an electronic filing cabinet, editing and printing output documents, supporting the exchange of information.

10.6 Monitoring of the technical condition is carried out as part of the diagnostic equipment within the time limits established in the medical institution.

11. Information about disposal in the territory of the Russian Federation.

11.1 Disposal of the product must be carried out in accordance with the standards in force on the territory of the Russian Federation - SanPiN 2.1.3684-21.

11.2 According to SanPiN 2.1.3684-21 the product belongs to class A - epidemiological safe waste.

11.3 The product must be disposed of if:

- end of service life;

- confirmation of facts and circumstances that endanger the life and health of health workers and indicate that the prescribed appointment has not been fulfilled.

11.4 Recycling is also subject to all packaging, including transport. Paper, wood, polyethylene and plastic must be disposed of separately in accordance with the rules for the collection, accounting and disposal of packaging.

12. Manufacturer's warranties

12.1 The manufacturer guarantees that the quality of the EEG systems meets the requirements of these technical specifications, subject to the rules of operation, transportation and storage.

12.2 Warranty period - 12 months from the date of commissioning.

Guaranteed shelf life - 6 months from the date of manufacture.

12.3 If there are justified claims to the quality of EEG systems, the manufacturer ensures their replacement during the warranty period of operation.

Claims for products with external damage will not be accepted.

12.4 Post-warranty repair is carried out for a separate fee.

12.5. The total service life / operation of the system is 10 years

The official representative of the manufacturer , as well as data for sending complaints and appeals :
Limited Liability Company "Infomed-Neuro", legal. address: 125445 , Moscow , Leningradskoe highway, 69, building 1, room 1 , floor 1, room 50 . phone 8 (495) 645-47-00 . kitaveva@im-neuro.ru .

Appendix A

WARRANTY CARD

Product category _____

Manufacturer _____

Model _____ Stamp _____

Serial number/art. _____ about sale _____

Date of sale _____

Organization of the seller _____

Seller's signature _____

The product has been tested and has no damage.

I have read and agree to the terms of use.

Buyer's signature _____

The warranty period is 2 years from the date of sale from the manufacturer's warehouse.

For warranty questions please contact:

Limited Liability Company "Infomed-Neuro", legal address: 125445 , Moscow , Leningradskoe highway, 69, building 1, room 1 , floor 1, room 50 . phone 8 (495) 645-47-00 . kitayeva@im-neuro.ru.

ATTENTION! Make sure that the warranty card is filled out correctly. The serial number of the product / article must exactly match the number in the coupon. All columns of the coupon (with the exception of the columns for the guarantee center) must be filled in.

In case of inaccurate or incomplete filling of the coupon, the warranty for the product is considered lost.

The manufacturer's warranty is carried out in accordance with the Civil Code of the Russian Federation on a general basis.

The product must not bear traces of mechanical or thermal effects, as well as traces of exposure to various chemical agents. The product may only be operated in accordance with the instruction manual for the product.

WARRANTY STATEMENT

DATE	DESCRIPTION	WARRANTY CENTER

Annex B

Specifications and specifications of accessories:

EEG electrodes are divided into two types:

-cranial:

Monopolar needle electrode

Electrode (plate) EL CFM for 2-channel EEG

Disc electrode with cable

Cabled electrodes for LTM

Pediatric disc electrode with cable

- ear

Ear electrode (clip) EL TPEAR (used with the EL TPLINK cable to combine two ear electrodes EL TPEAR)

Disposable electrodes EL PUD

Collecting surface diameter:

- cranial electrode no more than 10 mm

- ear electrode - 10 mm.

Electrical strength of insulation, V, not less than 30

Insulation resistance, Ohm, not less than 10^9

Electrode potential difference, mV, no more than 100

Drift of electrode potential difference (voltage drift), μV , not more than 25 for EEG electrodes

Electrode impedance, Ohm, no more than $5 \cdot 10^3$

Ready time, min, no more than 10

Noise voltage, not more than 20 μV for EEG electrodes

A non-removable electrode wire (cable) of reusable electrodes, stretched with a force of 1 N, must withstand at least 500 bends by $\pm 90^\circ$ at the point where the wire exits the electrode.

Disposable ECG electrode ECG EL, Cable for ECG electrode

Continuous contact time, h, not less than 0.5 for ECG electrodes

Noise voltage, not more than 30 μV for ECG electrodes

ECG electrode cables have vibration strength after exposure to vibration loads with a frequency of 30 Hz and an amplitude of 0.3 mm.

Drift of electrode potential difference (voltage drift), μV , not more than 250 for ECG electrodes

1. Electrode cap:

- MRINEO19 X

- MRINEO31 X

- MRINEO61 X

CAP MRINEO19X

19 electrodes + ground + pad

2 bipolar ECG/EMG inputs, 4 m cable, male plug, 37 pins

31 electrode + ground + pad

	1 bipolar ECG/EMG input, 4 m cable, male plug, 37 pins		
CAP MRINEO61X	61 electrode + ground + pad 1 bipolar ECG/EMG input, 4 m cable. male plug, 37 pins		
SD junction box should be used MRI - SD MRI 64			
FOR ADULTS	FOR ADULTS	FOR ADULTS	FOR ADULTS
x-7 - 58-62 cm	x=6 - 54-58 cm	x-5 - 50-54 cm	x \u003d 4 - 46-50 cm
big	average	small	very little
•	•	•	•
CAP MRINEO19 7	CAP MRINEO19 6	CAP MRINEO19 5	CAP MRINEO19 4
CAP MRINEO31 7	CAP MRINEO31 6	CAP MRINEO31 5	CAP MRINEO31 4
CAP MRINEO61 7	CAP MRINEO61 6	CAP MRINEO61 5	X

There is a reference (REF) and ground (GND) (spare) electrode.

Electrode impedance test: 1-2 min

2. Cable for EL CFM electrodes:

- CB CFM050

- CB CFM150

CB CFM050	Connecting cable for EL CFM electrodes, length 50 cm, with touch-proof plugs for each electrode (C3, C4, REF , GND)	
CB CFM0 1 50	Connecting cable for EL CFM electrodes, length 1.5 m, with touch-proof plugs for each electrode (C3, C4, REF , GND)	

3. Electrodes with cable for LTM :

- CB EL TPCO 19

- CB EL TPCO 31

BRAIDED CUP ELECTRODES (AG) WITH LABELED TOUCH-PROOF CONNECTORS 150 CM LONG + GROUND AND PAD FOR LONG-TERM MONITORING		
	CB EL TPCO 19	Cable with 19 electrodes + ground + pad
	CB EL TPCO 31	Cable with 31 electrodes + ground + pad

4. Electrode (plate) EL CFM for 2-channel EEG:

	EL CFM	Box with single electrodes for 2-channel EEG (CFM) + grounding + pad
---	--------	--

5. Disc electrode with cable:

- EL TPCO
- EL TPCO 6
- EL TPCO G06
- EL TPCO G10

	Cup electrodes (Ag), cable length - 90 cm	
	EL TPCO	Cup electrode 10 mm long with touch-proof connector
	EL TPCO 6	Cup electrode 6 mm long with touch-proof connector
	Cup electrodes with gold-plated contacts, cable length - 90 cm	
	EL TPCO G10	Cup electrode 10 mm long with touch-proof connector
	EL TPCO G06	Cup electrode 6 mm long with touch-proof connector

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WITH SILVER-PLATED	EL TPCO S 100	Cast silver cup, 10 mm, cable length - 100 cm	
	EL TPCO S 150	Cast silver cup, 10 mm, cable length - 150 cm	
	EL TPCO S 200	Cast silver cup, 10 mm, cable length - 200 cm	
	EL TPCO S 250	Cast silver cup, 10 mm, cable length - 250 cm	

WITH GOLD-PLATED	EL TPCO SG 100	Cast gold cup, 10 mm, cable length - 100 cm
	EL TPCO SG 150	Cast gold cup, 10 mm, cable length - 150 cm
	EL TPCO SG 200	Cast gold cup, 10 mm, cable length - 200 cm
	EL TPCO SG 250	Cast gold cup, 10 mm, cable length - 250 cm

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WITH SILVER-PLATED	EL TPCO SLT 100	Cast silver cup, 10 mm, cable length - 100 cm
	EL TPCO SLT 150	Cast silver cup, 10 mm, cable length - 150 cm
	EL TPCO SLT 200	Cast silver cup, 10 mm, cable length - 200 cm
	EL TPCO SLT 250	Cast silver cup, 10 mm, cable length - 250 cm
WITH GOLD-PLATED	EL TPCO SLTG 100	Cast gold cup, 10 mm, cable length - 100 cm
	EL TPCO SLTG 150	Cast gold cup, 10 mm, cable length - 150 cm
	EL TPCO SLTG 200	Cast gold cup, 10 mm, cable length - 200 cm
	EL TPCO SLTG 250	Cast gold cup, 10 mm, cable length - 250 cm

6. Pediatric disc electrode with cable:

- EL TPCO S6 100
- EL TPCO S6 150
- EL TPCO S6 200
- EL TPCO S6 250


- EL TPCO SG6 150
- EL TPCO SG6 200
- EL TPCO SG6 250

WITH SILVER-PLATED	EL TPCO S6 100	Cast silver cup, 6 mm, cable length - 100 cm	
	EL TPCO S6 150	Cast silver cup, 6 mm, cable length - 150 cm	
	EL TPCO S6 200	Cast silver cup, 6 mm, cable length - 200 cm	
	EL TPCO S6 250	Cast silver cup, 6 mm, cable length - 250 cm	
WITH GOLD-PLATED	EL TPCO SG6 100	Cast gold cup, 6 mm, cable length - 100 cm	
	EL TPCO SG6 150	Cast gold cup, 6 mm, cable length - 150 cm	
	EL TPCO SG6 200	Cast gold cup, 6 mm, cable length - 200 cm	
	EL TPCO SG6 250	Cast gold cup, 6 mm, cable length - 250 cm	

7. Ear electrode (clip) EL TPEAR with cable:

EL TPEAR	Pair of clip-on ear cup electrodes (AgCl) with 120 cm cable and touch-proof connector	
----------	---	---

8. EL TPLINK cable for connecting two EL TPEAR ear electrodes:

EL TPLINK	Connecting cable for two ear electrodes with contact pad (" G 2") and touch-proof connector at both ends	
-----------	--	---

9. Connecting cable for LTM / LTM 256:

	LTM SHORT	Adapter cable with touch-proof connector for SD junction box LTM , length - 40 cm, green (ground), black (for transmission of the reference signal) cable for shorting the input
	LTM SHORT 256	Adapter cable with touch-proof connector for SD LTM EXPRESS with 256-channel system, length - 40 cm, between each connector green (ground), black (for transmission of the reference signal) cable for shorting the input


10. CB PUD150 cable for EL PUD electrodes

11. Disposable electrodes EL PUD

	EL PUD	Box with 10 disposable electrodes for recording data from the external genitalia. Connecting cable required
	CB-PUD150	Connection cable for vulvar electrodes with 2 touch-protected stimulation connectors, with pentapolar electrode with DIN connector for identification

12. Monopolar needle electrode:

- EMG _ MBOT 2530
- EMG _ MBOT 2540
- EMG _ MBOT 3745
- EMG MBOT5050
- EMG MBOT7570

EMG MBOT 2530	length - 25 mm, diameter - 0.30 mm	●	
EMG MBOT 2540	length - 25 mm, diameter - 0.40 mm	●	
EMG MBOT 3745	length - 37 mm, diameter - 0.45 mm	●	
EMG MBOT5050	length - 50 mm, diameter - 0.50 mm	○	
EMG MBOT7570	length - 75 mm, diameter - 0.30 mm	●	

14. Snoring sensor:

- EPM 1200
- EPM 1250

EPM 1200	Piezo-electric snore sensor (sensitive to skin vibration, insensitive to background noise) with touch-proof connector	
EPM 1250	Dynamic snoring sensor (signal that captures sound intensity, registration of respiratory rate and snoring) with touch-proof connector	

15. Respiration sensor kit EPM 1312, includes sensor (EPM 1387) and 4 tapes in different sizes (EPM 1340 + EPM 1341 + EPM 1342 + EPM 1343):

EPM 1340	Respiratory monitor strap kit including transducer (EPM s 1387) with 4 different size straps (EPM s 1340+ EPM s 1341+ EPM s 1342+ EPM s 1343)
----------	---

COMPONENTS

	EPMs 1387 _	Piezocrystalline breath sensor with touch-proof connector
---	-------------	---

EPMs 1340 _	Extra large elastic band (190x3.8 cm)	
EPMs 1341 _	Large elastic band (135x3.8 cm)	
EPMs 1342 _	Small elastic band (90x2.5 cm)	
EPMs 1343 _	Elastic band pediatric (45x2.5 cm)	

16. Pediatric breathing sensor EPM 1390

17. Pediatric airflow sensor EPM 1400

	EPM 1400	Piezoelectric breath control sensor, pediatric, reusable single-channel manifold with touch-proof connector
	EPM 1390	Miniature respiratory sensor, pediatric, with strap for fastening and connector with touch protection

18. Set of EPM 1439FM/ GM / HM disposable breath sensors , includes: EPM 1415 cable and 10 sensors adult 1426 EPMs / medium size 1427 EPMs / infant size 1428 EPMs)

EPM 1439FM	Set of disposable breath sensors including EPMs 1415 cable and 10 adult transducers
EPM 1439GM	Set of disposable breath sensors including EPMs 1415 cable and 10 transducers (medium size)
EPM 1439 HM	Set of disposable breath sensors including EPMs 1415 cable and 10 pediatric transducers

COMPONENTS

	EPMs 1415	Air flow sensor cable, single channel, with touch-proof connector
	EPMs 1426	Connector for 20 sensors (for adults)
	EPMs 1427	Connector for 20 sensors (medium size)
	EPMs 1428	Connector for 20 sensors (pediatric)

19. Air flow sensor EPM 1450

EPM 1450	Single-channel thermistor with airflow sensor with converter (EPMS 1450- S) and control panel (EPMS 1450-B)	
COMPONENTS		
EPMs 1450-B	EPC 1450 airflow sensor interface with touch-proof connector	

EPMs 1450- S	EPM 1450 airflow sensor thermistor converter	
--------------	--	---

20. Set of air flow and pressure sensors EPM 1475

EPM 1475	Thermistor with airflow and pressure sensor, combined sensor (EPMs 1475) and ThermoCan cannula set (EPMs 15807)	
COMPONENTS		
EPMs 1475	Thermistor with combined airflow and pressure sensor with touch-proof connector used for ThermoCan cannula set (EPMs 15802 and EPMs 15805)	
EPMs 15802	ThermoCan cannula set (5 pcs.), 60 cm (2 positions)	
EPMs 15807	Cannula set ThermoCan (5 pcs.), 210 cm (7 positions)	

21. Air pressure sensor and snoring EPM 14830

EPM 14830	Piezoelectric pressure sensor for airflow and snoring with transducer (EPMS 14830) and adult CO2 monitoring cannula set (EPMS 14805)	
COMPONENTS		
EPMs 1475	Piezo pressure transducer for airflow and snoring EPMs with touch-proof connector	
EPMs 15802	Set of cannulas for CO2 monitoring (5 pcs.) for adults (nasal and oral)	
EPMs 15803	Set of CO2 monitoring cannulas (5 pcs.) for teenagers (nasal)	
EPMs 15804	CO2 monitoring cannula set (5 pcs.) for pediatric (nasal)	
EPMs 15805	Set of cannulas for CO2 monitoring (5 pcs.) for adults (nasal)	

22. Body position sensor EPM 1563

EPM 1563	DC set of body position sensors including sensors (EP M s 1563) and elastic band (EP M s 1550)	
COMPONENTS		
EP M s 1563	Body position sensor DC converter EPM 1563 with touch-proof connector	
EP M s 1550	EPM 1520 Body Sensor Elastic Band	

23. EPM 1715 motion sensor kit (including piezoelectric transducer (EPM 1770) and pairs of ankle straps (EPMS 1730))

EPM 1715	Motion sensor kit including piezoelectric sensors (EP M s 1770) with touch-proof connector and a pair of elastic bands with pads (EP M s 1550)	
COMPONENTS		
EP M s 1770	Piezoelectric motion sensor EPM 1715 with touch-proof connector	
EP M s 1730	Pair of padded elastic bands for use with EPM 1715 pressure transducers	

24. Chest force sensor kit, including interface cable and 4 chest straps:

- EPM 9102 A for adults
- EPM 9102 C pediatric

25. Set of abdominal force sensor, including interface cable and 4 abdominal bands:

- EPM 9101 A for adults
- EPM 9101 C pediatric

ERM 9101 A	4 adult inductive abdominal sensor strips (EPMs 9002 - L 90, EPMs 9002- L 120), sensor cable (EPMs 9101)	
EPM 9101 C	4 strips of inductive abdominal sensors for pediatrics (EPM s 9002- L 40, EPMs 9002- L 60), sensor cable (EPMs 9101)	
ERM 9102 A	4 strips of inductive chest breathing sensors for adults (EPM s 9002- L 90, EPMs 9002- L 120), sensor cable (EPMs 9101)	
EPM 9102 C	4 strips of inductive sensors for pediatric chest breathing (EPM s 9002- L 40, EPMs 9002- L 60), sensor cable (EPMs 9102)	

COMPONENTS

CABLE + INTERFACES		
	EPMs 9101	Interface cable with touch-proof connectors - for inductive abdominal sensors
	EPMs 9102	Interface cable with touch-proof connectors - for inductive chest sensors

26. A set of abdominal and thoracic force sensors, including abdominal and thoracic interface cables, and 4 tapes:

- EPM 914x A for adults
- EPM 914x C pediatric
- EPM 915x A for adults
- EPM 915x C pediatric

ERM 914xA	4 inductive strips for adults (EPMs 9002 - L 90, EPMs 9002- L 120), cable 45 cm (EPMs 914x)	
ERM 914xS	4 inductive pediatric bands (EPM s 9002- L 40, EPMs 9002- L 60), 45 cm cable (EPMs 914x)	
ERM 915xA	4 inductive strips for adults (EPMs 9002 - L 90, EPMs 9002- L 120), cable 210 cm (EPMs 915x)	
ERM 915xS	4 inductive pediatric bands (EPM s 9002- L 40, EPMs 9002- L 60), cable 210 cm (EPMs 915x)	

COMPONENTS

	A SET OF INDUCTIVE TAPES FOR THE MORPHEUS LTM EXPRESS EEG SYSTEM AND / OR A SET FOR THE SD LTM EXPRESS EEG SYSTEM, INCLUDING 2 CABLES (YELLOW FOR THE ABDOMINAL SENSOR, BLUE FOR THE CHEST BREATHING SENSORS) AND 4 DIFFERENT TAPES	
	ERM 914x	Cable length 45 cm
	ERM 915x	Cable length 210 cm
PACK OF 2 INDUCTIVE STRIPS		
	EPM s 9002- L 40	Small - 40 cm
	EPM s 9002- L 60	Average 60 cm
	EPM s 9002- L 90	Large - 90 cm
	EPM s 9002- L 120	Very large - 120 s
	EPM s 9002- L 133	Size XX L - 133 cm
	EPM s 9010	Pack of 1 length-adjustable inductive tape

27. Tape for breath sensor:

- EPMs 1340
- EPMs 1341
- EPMs 1342
- EPMs 1343

EPM s 1340	Extra large elastic band (190 x 3.8 cm)	
EPM s 1341	Large elastic band (135x3.8 cm)	
EPM s 1342	Small elastic band (90x2.5 cm)	
EPM s 1343	Elastic band pediatric (45x2.5 cm)	

28. Piezocrystalline breath sensor EPMs 1387

	EPMs 1387	Piezocrystalline breath sensor with touch-proof connector
---	--------------	---

29. Cable for airflow sensor single channel EPMs 1415

	EPMs 1415	Air flow sensor cable, single channel, with touch-proof connector
---	-----------	---

30. Thermal breath sensor:

- EPMs 1426 adult size
- EPMs 1427 medium size
- EPMs 1428 pediatric

	EPMs 1426	Connector for 20 sensors (for adults)
	EPMs 1427	Connector for 20 sensors (medium size)
	EPMs 1428	Connector for 20 sensors (pediatric)

31. EPMs 1450-B interface for EPM 1450 airflow sensor
S transducer for EPM 1450 airflow sensor

EPMs 1450-B	EPC 1450 airflow sensor interface with touch-proof connector	
EPMs 1450- S	EPM 1450 airflow sensor thermistor converter	

33. Combined thermistor pressure and air flow sensor EPMs 1475

EPMs 1475	Thermistor with combined airflow and pressure sensor with touch-proof connector used for ThermoCan cannula set (EPMs 15802 and EPMs 15807)	
EPMs 15802	Cannula set ThermoCan (5 pcs.), 60 cm (2 positions)	
EPMs 15807	Cannula set ThermoCan (5 pcs.), 210 cm (7 positions)	

CO 2 cannula :

- EPMs 14802, nasal and oral, male adult
- EPMs 14803, male for newborns
- EPMs 14804, male pediatric
- EPMs 14805, male adult

35. Piezoelectric pressure transducer for air flow and snoring EPMs 14830

EPMs 14830	EPMs 14830 Airflow and Snoring Piezo Sensor Transducer with Touch-Proof Connector	
EPMs 14802	Set of cannulas for CO2 monitoring (5 pcs.) for adults (nasal and oral)	
EPMs 14803	Set of CO2 monitoring cannulas (5 pcs.) for teenagers (nasal)	
EPMs 14804	CO2 monitoring cannula set (5 pcs.) for pediatric (nasal)	
EPMs 15805	Set of cannulas for CO2 monitoring (5 pcs.) for adults (nasal)	

36. Elastic tape EPMs 1550 for body position sensor EPM 1520

37. EPMs 1563 transducer for EPM 1563 body position sensor

EP M s 1563	Body position sensor DC converter EPM 1563 with touch-proof connector	
EP M s 1550	EPM 1520 Body Sensor Elastic Band	

38. ThermoCan cannula:

- EPMs 15802
- EPMs 15807

EPMs 15802	ThermoCan cannula set (5 pcs.), 60 cm (2 positions)	
EPMs 15807	Cannula set ThermoCan (5 pcs.), 210 cm (7 positions)	

39. Pair of EPMs 1730 ankle straps

40. Piezoelectric transducer EPMs 1770 for motion sensor EPM 1715

ERM 1715	Motion sensor set including piezoelectric sensors (EP M s 1770) with touch-proof connector and a pair of elastic bands with pads (EP M s 1730)	
COMPONENTS		
EP M s 1770	Piezoelectric motion sensor EPM 1715 with touch-proof connector	
EP M s 1730	Pair of padded elastic bands for use with EPM 1715 pressure transducers	

41. A Pair of Induction Tapes:

- 9002- L 120 EPMs
- 9002- L 133 EPMs
- 9002-L40 EPMs
- 9002-L60 EPMs
- 9002-L90 EPMs

	EPM s 9002- L 40	Small - 40 cm
	EPM s 9002- L 60	Average 60 cm
	EPM s 9002- L 90	Large - 90 cm
	EPM s 9002- L 120	Very large - 120 s
	EPM s 9002- L 133	Size XX L - 133 cm

42. Adjustable induction band EPMs 9010

	EPM s 9010	Pack of 1 length-adjustable inductive tape
---	------------	--

- EPMs 9101 for the abdomen
- EPMs 9102 for the thoracic region

	EPMs 9101	Inductive interface cable with touch-proof connectors - for inductive abdominal sensors
	EPMs 9102	Inductive interface cable with touch-proof connectors - for inductive chest breathing sensors

44. Induction interface cable with connector for MORPHEUS:

- EPMs 914x
- EPMs 915x

EPMs 914x	Cable length 45 cm
EPMs 915x	Cable length 210 cm

45. Disposable ECG electrode ECG EL

ECG EL	Monopolar self-adhesive ECG electrode with 26 mm disposable clip, pack of 30, min. 5 pieces.	
--------	--	--

Overall diameter 26 mm, collector surface diameter 24 mm

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
Annex A1


Drawings and diagrams

Correct connection of various devices


grounding to ground; isolation transformer


DATA COLLECTION SYSTEM; POINT-TO-POINT CONNECTION (LAN CABLE); VIEW STATION;


grounding to ground; isolation transformer


PATIENT ZONE; Micromed system; BQ USB or BQ USB Express; Overlay part of Micromed; Network separator; Micromed system; BQ PCI EXPR LTM HALF/FULL, BQ USB or BQ USB Express; Overlay part of Micromed; isolation transformer; Micromed system; ETHERNET; Micromed system; EXPRESS TERMINAL BOX; Overlay part of Micromed; BQNET POE 2/M; Overlay part of Micromed; Electrical signal channel; network connection; Processing and storage of digital signals; Connection; Signal processing and digitization


Scheme.- cable connection via interface BQ PCI EXPR LTM HALF / FULL in the patient's room


Scheme - cable connection via interface BQ PCI EXPR LTM HALF / FULL and PC outside the patient room



Wiring diagram for long distances


Scheme - Configuration 128 – channel system with two amplifiers SD LTM 64 EXPRESS


Scheme – the location of the amplifier inputs JB LTM 32 EXPRESS

Scheme - the location of the inputs of the amplifier JB LTM 64 EXPRESS

amplifier block diagram

Rice. – cable connection with BQ USB EXPRESS





Diagram - cable connection via BQ interface PCI EXP LTM HALF / FULL in the room with the patient

Schematic - 128-channel system configuration with two SD amplifiers LTM 64 EXPRESS


JB input layout LTM 32 EXPRESS


JB input layout LTM 64 EXPRESS

Figure 2 - Cable connection with BQ USB


Scheme - cable connection with BQ PCI interface


Scheme. Main parts of the device


Figure 2 (above) shows the location of the main components of the device, and below - the location of the connectors and a description of their functions. For further information on connector specifications, see section " Mechanical Specifications " on page 30.


Rice. 3. Diagram of input connectors and controls


Block diagram BQ EXPRESS USB connection


Block diagram BQ PCI EXPRESS connection next to the patient


Block diagram Connection BQ PCI EXPRESS and PC outside the patient area

Hbc/ 8/ BOX TERMINAL EXPRESS


Scheme. Stretch Connection Diagram


Wiring Diagram of Gain Tests

BOX TERMINAL EXPRESS

KVM wiring diagram


- 1 Source (computer, processor)
- 2 KVM module
- 3 Connection cable
- 4 KVM Extender
- 5 Console (monitor, keyboard, mouse)
- 6 Second monitor (optional)
- 7 USB devices


KVM chassis (images). It is important to understand that the difference in
2 names directly depends on the number of slots)

KVM Chassis KVM BODY 2


1. Slot for KVM modules
2. Connecting to a 5V DC Power Supply


KVM Chassis KVM BODY 6 R

- 1,2,3 Slot for KVM modules
- 4 connection to a standard power supply
- 5 Connect to a 5V DC power supply


KVM Chassis KVM BODY 4 R

- 1,2 Slot for KVM modules
- 3 connection to a standard power supply
- 4 Connect to 5V DC power supply

KVM Chassis KVM BODY 6

- 1,2,3 Slot for KVM modules
- 4 Connect to 5V DC power supply

KVM modules (images and description)

Rice. 1


Fig 2


Fig.3

Fig.1 description (KVM - DVI module TX AUD dual, DVI TX AUD DCAT dual, DVI RX AUD)

No.	Description	Function
1	grounding	the ground wire (used to ground the unit) is connected.
2	KVM port	KVM cable is connected that connects the module to your computer. Only KVM cables can be used
3	Remote LED	There are two LEDs on the local unit to indicate the operating status of the local and remote units.
4	Local LED	
	Working mode Choice	Pressing this button toggles between local and remote modes of operation.

Rice. 2 description (DVI KVM module RX , DVI TX DCAT ,

No.	Description	Function
1	grounding	the ground wire (used to ground the unit) is connected.
2	Line LED	These LEDs show the operating status of local and remote devices.
3	Power indicator	

Rice. 3 description (KVM module KIT KVM, KVM module KIT KVM DUAL)

No.	Description	Function
-----	-------------	----------

1	Power connector	Connect the cable from the DC power adapter to this connector. *
2	Mouse port	Console mice connect to these ports.
3	Keyboard port	Console keyboards connect to these ports.
4	Link Port	The category 5 cable connecting the remote and local units is
5	Video port	Console monitors connect to these ports.

*Note: The module can be powered by a computer - external power is only required if the power from the local computer/computers in the KVM installation is insufficient.

Scheme of connecting KVM to PC

1 Connect the cables from the local console devices (keyboard, monitor, mouse) to their ports on the back of the local device (CE250 **AL**).

2 Connect the SPHD end of the KVM cable supplied with this unit to the KVM port on the front panel of the CE250. **AL** .

Note. The shape and function of the connectors on the cable and the KVM switch have been changed so that only KVM cables designed to work with this switch can be used.

3 Connect the connectors on the other end of the cable to the appropriate ports on your computer (or console section of the KVM switch if you are using one). Each connector is marked with a corresponding icon to identify itself.

4 Connect one end of a category 5 or higher STP cable to the CE250. **REMOTE I/O AL** port. Connect the other end of the cable to the CE250. **AR Remote I/O** port.

Note. The STP cable is not included in the package. Separate purchase required. Cable length can be up to 150 m (500 ft).


5 Connect the cables from the remote console devices (keyboard, monitor, and mouse) to their ports on the CE250 console side. **AR** .

6 Connect the power adapter (supplied) to AC power; connect the power cable of the adapter to the remote unit (CE250 **AR**) Power socket.


7 If you choose to use the power adapter with the CE250, **PLEASE** connect the power adapter to an AC power source; Connect the power cable of the adapter to the CE250 **AL** Power Connector.

Note: The use of a power adapter with the CE250 **AL** is optional. See Power connector note on page 7 for details.


115


BQ CONTACT contact closure control module (optional)

CAMERA IP COMPACT video camera (optional)

- 1 focus lock ring
- 2 Focus ring
- 3 IR LED
- 4 PIR sensor
- 5 speaker
- 6 Status indicator
- 7 Control button
- 8 Not used
- 9 PoE network connector
- 10 SD card slot
- 11 Power connector


CAMERA IP HD OUTDOOR video camera (optional)


Camera Terminal block

Terminal Block Diagram

CAMERA IP HD PRO video camera (optional)

1 SD card slot (microSD)

2 Control button

3 Status LED

4 RJ45 connector


5 Audio input

6 Audio output


7 I/O connector

8 Power connector

Camcorder CAMERA IP HD PLUS





- 1 power connector
- 2 3.5mm unbalanced audio input
- 3 3.5 mm unbalanced audio output
- 4 Network connector
- 5 Microphone power indicator (XLR)
- 6 XLR. balanced audio input (left)
- 7 XLR. balanced audio input (right)
- 8 BNC connector (SDI)
- 9 HDMI connector

- 10 SD card slot (SDHC)
- 11 Status LED
- 12 I/O connectors
- 13 Control button
- 14 indicator lights (3x)

NET ISOLAT Network Isolator for 10/100 Mbps Networks

Left

Right

Figure 11. Schematic diagram of the EN-1005+.

Annex A2

Data on electromagnetic compatibility of system parts

BRAIN QUICK systems suitable for use in the specified electromagnetic environment. The user must ensure that the device is used in an electromagnetic environment that is compatible with the following:		
Radiation Tests	Correspondence	Electromagnetic environment
Conducted and radiated RF emissions CISPR 11	class B	BRAIN QUICK systems are suitable for use in all establishments except residential buildings and those directly connected to public low-voltage networks.
Conducted and radiated RF emissions CISPR 11	group 1	BRAIN QUICK systems use RF energy only for their own functioning. Therefore, RF emissions are very low and should not interfere with nearby electronic equipment.
Harmonic emissions IEC 61000-3-2	Compliant with class A	BRAIN QUICK systems are suitable for use in all establishments except residential buildings and those directly connected to public low-voltage networks.
Voltage fluctuations / flickering radiation IEC 61000-3-3	Corresponds	BRAIN QUICK systems are suitable for use in all establishments except residential buildings and those directly connected to public low-voltage networks.

Stability tests	Test level EN 60601-1-2	Compliance level	Electromagnetic environment
Electrostatic discharge (ESD) EN 61000-4-2	8 kV contact 2/4/8/15 kV air	Not applicable	Floors should be covered with wood, concrete or ceramic tiles. If the floors are covered with synthetic material, the relative humidity must be at least 30%.
Fast electrical transients/bursts IEC 61000-4-4	2 kV for power lines 1 kV for I/O lines > 3 m	EN 60601-1-2 Test level	Mains power quality should be that of typical commercial and hospital environments.
throws IEC 61000-4-5	0.5/1 kV differential mode 0.5/1/2 kV common mode	EN 60601-1-2 Test level	Mains power quality should be that of typical commercial and hospital environments.

Registrar **SD LTM 32/64 EXPRESS** designed to operate in special electromagnetic environments. Users **SD LTM 32/64 EXPRESS** must ensure that the device used in the electromagnetic environment complies with the following conditions:

Emission Tests	Correspondence	Electromagnetic environment
----------------	----------------	-----------------------------

Conducted and Emissive CISPR radiation eleven	RF group I	SD LTM 32/64 EXPRESS uses RF energy only for its work. Therefore, the RF emission is very weak and does not cause interference in the surrounding
	Class B	
Harmonic radiation IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker IEC 61000-3-3	Not applicable	SD LTM 32/64 EXPRESS can be used in all installations, except for domestic and / or those directly connected to the low voltage network.

Immunity Tests	Test level	Compatibility level	Electromagnetic environment I
Electrostatic Discharge (ESD) EN 61000-4-2	Contact 8 kV 2/4/8/15kV on air	Not applicable	Floors should be wood, concrete or ceramic. If the floors are covered with synthetic material, the relative humidity must be at least 30%.
Short-term instability in the power supply / impulses IEC 61000-4-4	2kV for power lines 1kV for input/output lines > 3 m	Not applied	The power supply must meet the power requirements of commercial buildings or hospitals.
Inrush current and voltage IEC 61000-4-5	0.5/1kV differential mode 0.5/1/2kV common	Not applicable	The power supply must meet the power requirements of commercial buildings or hospitals .
Voltage drops (dips), short-term outages and voltage fluctuations in power networks IEC 61000-4-11	0% Un for 0.5 cycle 0 % Un for 1 cycle, 70 % Un for 25 cycles 0 % Un For 5 sec	Not applicable	The power supply must meet the power requirements of commercial buildings or hospitals . If continuous operation is required in the event of power outages, it is recommended to use a UPS or battery.
Mains frequency (50/60 Hz) IEC 61000-4-8	thirty A/m	EN 60601-1-2 Test level	The frequency characteristics of the current must be suitable for the power supply requirements of commercial buildings or hospitals.

RF emissions EN 61000-4-3	3V/m 80 MHz up to 2.7GHz	EN 60601-1-2 Test level	Portable and mobile RF communications equipment is not must be near the recorder, including cables . The minimum distance is 30 cm.
Conducted RF emissions EN 61000-4-6	3B 150 kHz -80 MHz		

Registrar **MORPHEUS** designed to operate in special electromagnetic environments. Users **MORPHEUS** must ensure that the device used in the electromagnetic environment complies with the following conditions:

Emission Tests	Correspondence	Electromagnetic environment
Conducted and Emissive RF Radiation CISPR eleven	group I	MORPHEUS uses RF energy only for its work. Therefore, the RF emission is very weak and does not cause interference in the surrounding electronic equipment. MORPHEUS can be used in all installations that are directly connected to the domestic low voltage network.
	Class B	
Harmonic radiation IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker IEC 61000-3-3	Not applicable	

Immunity Tests	Test level EN 60601-1-2	Compatibility level	Electromagnetic environment ¹
Electrostatic Discharge (ESD) EN 61000-4-2	Contact 6 kV 8kV in air	Test level EN 60601-1-2	Floors should be wood, concrete or ceramic. If the floors are covered with synthetic material, the relative humidity must be at least 30%.
Short-term instability in the power supply / impulses IEC 61000-4-4	±2kV for power lines ±1kV for input/output lines	Test level EN 60601-1-2	The power supply must meet the power requirements of commercial buildings or hospitals.
Surge current and voltage IEC 61000- 4-5	±1kV differential mode ±2kV common	Not applicable	The power supply must meet the power requirements of commercial buildings or hospitals .

Voltage drops (dips), short-term outages and voltage fluctuations in power networks IEC 61000-4-11	<5% U_t² (>95% drop in U_t) for 0.5 cycle, 40% U_t (60% drop in U_t) for 5 cycles, 70% U_t (30% drop in U_t) for 2.5 cycles, <5% U_t (>95% drop in U_t) for 5 sec	Not applicable	The power supply must meet the power requirements of commercial buildings or hospitals . If continuous operation is required in the event of power outages, it is recommended to use a UPS or battery.
Mains frequency (50/60 Hz) IEC 61000-4-8	3 A/m	Test Level EN 60601-1-2	The frequency characteristics of the current must be suitable for the power supply requirements of commercial buildings or hospitals.

RF emissions EN 61000-4-3	Auxiliary equipment	EN test level 60601-1-2	Portable and mobile RF communications equipment is not should be next to the MORPHEUS registrar , including cables. The recommended distance is calculated based on the frequency of the transmitter. Recommended distance: $d = 1.2\sqrt{P}$ from 150kHz up to 80MHz ³ $d = 1.2\sqrt{P}$ from 80MHz to 800MHz ³ $d = 2.3\sqrt{P}$ from 800MHz up to 2.5GHz Where P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended distance in meters (m). Magnetic fluxes from RF transmitters determined from field studies ⁴ , should not exceed the permissible level in each frequency range. ⁵  Interference may occur in the vicinity of equipment marked with the following symbol:
Conducted RF* EN 61000-4-6	3V/m 80 MHz - 2.5GHz 3Vef [10Veff in range ISM] 150kHz -80MHz		

• not applicable

Comment 1: these recommendations cannot be applied in every situation. The propagation of electromagnetic radiation is influenced by the processes of absorption and reflection from structures, objects and people.

Comment 2: UT is the mains AC voltage before the test level is applied.

Comment 3: at frequencies 80MHz and 800MHz should apply the high frequency range equation.

Comment 4: magnetic flux from transmitters, e.g. base stations for radio (cellular/wireless) telephones and land mobile radios, amateur radios, broadcast radios AM and FM television is theoretically unpredictable. To assess the electromagnetic environment using RF transmitters, a site survey must be performed. If the measured magnetic flux at the location of **MORPHEUS** exceeds the permissible RF level, the correct operation of the recorder should be monitored **MORPHEUS**. In the event of a malfunction, additional measures should be taken, such as relocating the recorder **MORPHEUS**.

Comment 5: at frequencies above 150 kHz up to 80MHz magnetic flux should be lower 3 V/m.

These characteristics are provided under the following conditions:

MORPHEUS equipped with shorter electrode cables 1.5m.

The connection to the PC is made using a cable, as indicated in the relevant section, or via radio. The use of accessories, transmitters, and cables other than those recommended by the manufacturer will increase emissions and reduce the immunity of **MORPHEUS**.

Electromagnetic environment

The SD Plus reception amplifier is suitable for use in the electromagnetic environment specified. The purchaser or user of SD Plus must ensure that the device is used in the electromagnetic environment specified below:

Radiation test	Correspondence	Electromagnetic environment
Emitted or conducted radio emission CISPR 11	Class B	SD device Plus Suitable for use in any establishment other than residential and/or directly connected to a low voltage power supply.
Emitted or conducted radio emission CISPR 11	Group I	This unit is SD Plus uses RF energy only for its internal functions. Therefore, RF emissions are negligible and are unlikely to interfere with nearby electronic equipment.
Harmonic radiation IEC 61000-3-2	Not applicable	SD Plus device Suitable for use in any establishment other than residential and/or directly connected to a low voltage power supply.
Voltage fluctuations / spike radiation IEC 61000-3-3	Not applicable	SD Plus device Suitable for use in any establishment other than residential and/or directly connected to a low voltage power supply.

Stability test	Test level EN 60601-1-2	Compliance Level	Electromagnetic environment requirements
Electrostatic discharge EN 61000-4-2	8 kV contact 2/4/8/15 kV air	Not applicable	The floor must be wood, cement or ceramic tiles. If the floor is covered with synthetic material, the relative humidity must be at least 30%
Radiated RF EN 61000-4-3	3 V/m 80 MHz to 2.7 GHz	Test level EN 60601-1-2	Portable and mobile radio communications equipment must be used at least 30 cm away from any part of the device, including cables
Induced RV EN 61000-4-6	3 V 150 kHz - 80 MHz		
Fast electrical transients or spikes EN 61000-4-4	2 kV for power lines 1 kV for input/output lines > 3 m	Not applicable	Mains power quality should be that of a typical commercial or hospital environment
overvoltage EN 61000-4-5	0.5/1 kV differential mode 0.5/1/2 kV common mode	Not applicable	Mains power quality should be that of a typical commercial or hospital environment

Voltage dips, momentary voltage interruptions and voltage fluctuations in power lines EN 61000-4-11	0% U_n for 0.5 cycle 0 % U_n for 1 cycle 70% U_n for 25 cycles 0 % U_n for 5 s	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If continuous use is required also during a power failure, it is recommended to power the system from an uninterruptible power supply or batteries
Magnetic field with mains frequency (50/60 Hz) EN61000-4-8	30 A/m	Test level EN 60601-1-2	Mains frequency magnetic field should be of a typical commercial or hospital location

Annex B
SCROLL

equipment, instruments and materials necessary to control the EEG system

Name	Main characteristics
Adjusting laboratory autotransformer	187...242 V, 8 A
Universal measuring device	0.06...1500 mA, 0.075...600 V, 1.5...600 V, 0.3...5000 kΩ, 0.5 μF
Ammeter	0.1-50 A, class 0.5
Voltmeter universal	0...2 A, 0...1000 V, ±0.5%
Load equivalents	-
Electronic frequency counter	0.1 Hz...120 MHz, $\pm 6 \times 10^{-7}$
Electronic oscilloscope	0-2 MHz, ±5%
ampervoltmeter	~0.05...750 V, 0.0002...7.5 A, class 1, 0, 45...300 Hz
Megohmmeter	0 to 100 MΩ (±1.5%)
microohmmeter	0...100 mΩ, ±2.5%
Breakdown setup	≈0...10 000 V, 0.5 kV · A, +4%
Two-hand stopwatch	0...60 s (0...30 min.), ±0.5 s
Climate chamber	from minus 70 °C to plus 90 °C ±2 °C; 10-100% at 10-60 °C, ±3%
Electric wall clock	
Resistor	1 kOhm, ±1%
Capacitor	0.15uF, ±5%
Vibrostand electrodynamic	5...60 Hz, ±0.5 Hz 0.1- 6 mm, ±15%
Impact stand	30...300 m/s ² , ±20% 1...16 m/s, 50...180 beats/min.
Sound level meter	20 Hz...100 kHz, A, B, C, D, lin. 24-170 dB, cl. 0
Ruler metal	0...1000 mm, ±1.0 mm
Scales	2.5...150 kg, division value 0.05 kg
RFI meter with network equivalent NNB101	0.1...30 MHz. from minus 30 to plus 132 dB

Note - It is allowed to use other control and measuring equipment and auxiliary equipment, provided that their technical characteristics are not worse than the above.

Annex B.1.

(reference)

Reference Documentation List

European standards

Art. No.	Title:	Device classification	Conformity assessment procedure
93/42/CEE	COUNCIL DIRECTIVE 93/42/ EEC of 14 June 1993 concerns medical devices	Class IIb	Application II
2011/65/EU	Directive 2011/65/ EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (modified)	--	--
1999/5/EC	Directive of the European Parliament and the EU Council of 9 March 1999 on radio equipment and telecommunications equipment and the mutual recognition of their conformity	--	--

All harmonized standards are listed below.

Agreed Standard	Name:
EN 60601-1-1:1990 + \A1+A11+A12: 1994 +\A2: 1997	Medical electrical equipment – Part 1: basic safety requirements
EN 60601-2-26:2003	Medical electrical equipment. Part 2-24. Particular safety requirements for electroencephalographs
EN 60601-1-2:2001	Medical electrical equipment. Part 1-2. General safety requirements. Additional standard. Electromagnetic compatibility. requirements and tests.
EN 60601-1-6:2010	Medical electrical equipment - part 1: basic safety requirements and important features – Substandard : applicability
EN 62366 : 2008	The equipment is medical. Use of technologies as applicable to medical equipment
EN ISO 14971:2012	Medical equipment – Application of risk management to medical devices
EN 62304 :2006-07	Medical equipment software. Software Product Life Cycle
ETSI EN 300 328 V1.9.1	Electromagnetic Compatibility and Radio Controlled Models (ERM): broadband data transmission systems; equipment transmitting data at a frequency of 2.4 GHz, modulation methods in the ISM radio frequency band and broadband band; agreed EN and containing the essential requirements of article 3.2 of the R & TTE directive

* These products are designed and manufactured by Micromed Quality System according to EN standards ISO 9001 and EN ISO 13485 (Certificate No. 9120, MIC 3 and 9124, MIC 4 issued by CSQ) and approved by Notified Body IMQ number 0051. Any change made to the product is carried out in strict accordance with the requirements of the Quality System and in accordance with the previously noted requirements. It is planned to take the necessary corrective actions based on the experience of using devices (equipment) in the operational cycle. All quality records related to products are kept for at least 5 years after the last product was manufactured.

Annex D

(mandatory)

Technical documentation for the components of the system

(see attached files and information)

VMIXER PRO Composite Switcher with 2 inputs / VMIXER PRO Professional Switcher with 5 HDMI inputs (optional) (C 2-6104 A Multi-window processor - article in the manufacturer's QMS. the article is indicated in the operational documentation)

Computer input

- DVI-D (on DVI-I connector) supporting up to 162 MHz (1920x1200 @ 60 Hz).
- HDCP support.
- EDID 1.3 data on DVI-D connector to convey available resolutions.
- Analog RGB/YPbPr supporting RGBHV, RGBS, RGsB, YPbPr, 75 ohm auto termination
- Digital Sync (in RGBHV or RGBS mode): TTL level, 10K interrupt, positive or negative
- Analog Sync (in RGsB, YPbPr, YUV modes): 0.3 V negative.
- Analog RGB Level Range: Approx. 0.5–2.0 Vp-p.
- Scan speed detection: automatic
- Analog PC Resolution: Any up to 2048x2048
- HDTV analog resolutions: any up to 1080p
- Analog maximum horizontal frequency: 150 kHz

Input type	5x DVI-U / HDMI via DVI-I connector
DVI-D Digital Format:	TMDS (Transition Minimized Differential Signaling)
Maximum HD Resolution	1080p at 60Hz
Maximum PC Resolution	1920x1200 @ 60Hz
DVI-A analog format:	Composite, S-Video, YPbPr, RGBHV
TV standards	NTSC, PAL, PAL-M, PAL-N, SECAM
Maximum HDTV Resolution	1080p at 60Hz
Maximum PC Resolution	2048x2048
Maximum Horizontal Scan Speed	150 kHz
Signal Format (via DVI-I)	RGBHV, RGBS, RGsB, YPbPr, YC, CV
Synchronize	TTL level, 10KQ, positive or negative
Termination	75Q
RGB level range	0.5-2.0 Vpk-pk
Scan speed detection	Automatically
DVI-D digital format	1x DVI-I / HDMI via DVI connector

Maximum HD Resolution	TMD5
Maximum PC resolution	1080p at 60Hz
Analog format DVI-A Maximum HDTV resolution	1920x1200 @ 60Hz
Maximum PC resolution	1080p at 60Hz
Vertical Refresh Rate	2048x2048
Signal Format	Any up to 250 Hz
RGB level	RGBHV, RGBS, RGsB, YPbPr
General	0.7 Vpk-pk
HDMI audio support	Paste switch and follow window A
HDCP Compliant	YES
Size and position	Automatic via AutoSet or manually
Image size	User Presets
Image freeze	One video frame
Settings memory	non-volatile
Zoom range	Variable up to 10x magnification
Shrinkage range	Variable up to 10%
Image mirroring	Horizontally and/or Vertically
Horizontal Filtration	Digital
YPbPr deinterlacing	Adaptive motion + noise reduction
Film Mode (525i / 1080i)	3:2 descent detection
Transformation Technologies	Own - CORIO®2
Frame rate conversion	Adding/removing a frame
Color Resolution	24 bits (16.8 million colors)
Sample rate	162 MHz
Digital sampling	24 bit, 4:4:4 format
ROM	Flash, upgradeable via RS-232
Analog video setup	RGB level contrast, phase

KVM -|MODULES, EXTENSIONS, CHASSIS

Modules, extenders, chassis provide access to a computer system from a remote console (keyboard, monitor and mouse). For management and security purposes, since you can keep the system unit in a safe place while you put the console in a place that provides easy access.

Provides an ideal means of accessing a remote system unit using the compact Cat 5 cable commonly used in Ethernet networks. Cat 5 provides a neat, convenient and reliable data connection.

Key feature of the KVM module KIT : Built-in 8KV ESD protection and 2KV surge protection. It also features a dedicated ASIC chip for maximum reliability and compatibility. It detects the distance to the system and automatically adjusts the gain to compensate.

Power supply 5V DC.

They are stands for K VM - modules.

2-segment chassis

Block	145 x 147 x 44 mm (5.7 x 5.8 x 1.7 inches)
-------	--

4-segment chassis

Block	293 x 147 x 44 mm (11.5 x 5.8 x 1.7 inches)
-------	---

6-segment chassis

Block	442 x 147 x 44 mm (17.4 x 5.8 x 1.7 inches)
-------	---

Modules

All analog (VGA) or digital (DVI) signals conforming to the DVI-I Single Link standard can be transmitted. This includes, for example, monitor resolutions such as 1920x1200 at 60Hz, Full HD (1080p) or 2K HD (up to 2048x1152). The data rate is limited to 165 MPixel/s.

Power supply: AC 100-240V 50/60Hz

Modules:

- Analog audio / serial port: max. 300 mA
- Serial (RS422): max. 150 mA
- Balanced Sound: Max.500mA
- Digital audio: max. 300 mA
- USB HID: max. 300 mA

KVM module KIT

Function			Meaning
Console selection			1 x button
Connectors	Console ports	KB	1 x 6-pin Mini - DIN F (purple)
		video	1 x HDB-15 F
		Mouse	1 x 6-pin Mini - DIN F (green)

LEDs	KVM ports	1 x SPHD-15F
	Power	1 x DC Jack
	From unit to unit	1 x RJ-45
	Local	1 (green)
	Remote	1 green)
Emulation	Power	1 green)
	Online	1 green)
	Keyboard / Mouse	PS/2
video		1280 x 1024 @ 60Hz (150m) DDC ; DDC2 ; DDC 2B *
Power consumption		5.0 V DC; 1.25W (Max) or when using an adapter 100-240V 50/60Hz
Wednesday	Working temperature.	0-50°C
	Storage temperature.	-20-60°C
	Humidity	0-80% RH
Physical Characteristics	Frame	Metal;
	Weight	0.33 kg
	Dimensions (L x W x H)	12.50 x 8.15 x 2.50 cm

Analog audio:

Transmission format	Almost CD-quality digitized audio (16-bit, 38.4 kHz)
Signal level	Line level (maximum 5 volts pk-pk)
Input impedance	47 kOhm
CPU Unit Connections	2 x 3.5mm stereo jacks (audio in and audio out)
CON Unit Connections	2 x 3.5mm stereo jacks (audio in and audio out)

Analog Audio USB 2.0

Transmission format	Almost CD-quality digitized audio, 16-bit (8, 11.025, 16, 22.05, 32, 44.1, 48 kHz)
Signal level	Line level (maximum 5 volts pk-pk)
Input impedance	20 kOhm
CPU Unit Connections	1x USB type B female
CON Unit Connections	2 x 3.5mm stereo jacks (audio in and audio out)

Compatibility	AES/EBU, S/PDIF, EIAJ CP1201, IEC60958
Standards	Dolby Digital, DTS, PCM
Bit depth	24 bit
Sampling frequency	32 to 96 kHz
CPU block (inputs)	• Mini-XLR (AES/EBU; symmetrical, lockable) • Coaxial (S/PDIF; RCA, Cinch) • Optical (S/PDIF; TOSLINK)

CON Block (Outputs)	• Mini-XLR (AES/EBU; symmetrical, lockable) • Coaxial (S/PDIF; RCA, Cinch) • Optical (S/PDIF; TOSLINK)
----------------------------	--

Bit depth	24 bit
Sampling frequency	32 to 192 kHz
Input Level	Maximum. 6.4 dBu balanced (gain: 0 dB) Max. 0.4 dBu unbalanced (gain: 0 dB)
Output Level	dBu balanced dBu unbalanced
phantom power	48 VDC
pre-gain	10 - 65 dB
CPU Unit Connections	2 stereo jacks 6.35 mm (2 audio inputs)
CON Unit Connections	2x 6.35mm stereo jacks (2 audio outputs)

Connecting cable type **cat X**

Solid Cat Cable X AWG 24	S/UTP cable (Cat 5e) according to EIA/TIA-568-B. Four pairs of AWG24 wires. Connection according to EIA / TIA-568-B (1000BASE-T).
Cat Patch Cable X AWG 26 / 8	S/UTP cable (Cat 5e) according to EIA/TIA-568-B. Four pairs of wires AWG26 / 8. Connection according to EIA / TIA-568-B (1000BASE-T).

Maximum allowed cable length

Cat Installation Cable X AWG 24	140 m (400 ft)
Cat Patch Cable X AWG 26 / 8	70 m (200 ft)

Fiber. Connecting cable type.

cable type	Characteristics
Singlemode 2T:00	• PCI IV (ZN) H 2E9 (in-house patch cable) • IV (ZN) HH 2E9 (internal switching cable) • I / AD (ZN) H 4E9 (internal or external cable, resistant) • A / DQ (ZN) B2Y 4G9 (floor cable, rodent proof)
Multimode 50pm	• Two fibers 50pm • IV (ZN) H 2G50 (in-house patch cord) • I / AD (ZN) H 4G50 (internal or external cable, resistant)

Multimode 62.5 pm	- Two fibers 62.5 pm - IV (ZN) HH 2G62.5 (internal switching cable) - A / DO (ZN) B2Y 4G62.5 (floor cable, rodent proof)
-------------------	--

Maximum allowed cable length

cable type	Maximum allowed cable length
Singlemode 21:00	10,000 m (32,800 ft)
Singlemode 21:00 XV	5,000 m (16,400 ft)
Multimode 50pm (OM3)	1000 m (3280 ft)
Multimode 50pm	400 m (1300 ft)
Multimode 62.5 pm	200 m (650 ft)

Extension cords

Contact	Signal	Contact	Signal	Contact	Signal
1	TMDS data 2-	9	TMDS data 1-	17	TMDS Data 0-
2	TMDS 2+ data	10	TMDS data 1+	18	TMDS data 0+
3	Data TMDS 2 GND	eleven	TMDS Data 1 GND	19	TMDS Data 0 GND
4	NC	12	NC	20	NC
5	NC	13	NC	21	NC
6	DDC input (SCL)	14	+5VDC high resistance	22	TMDS Clock GND
7	DDC output (SDA)	15	GND	23	Hours TMDS+
8	Vertical Sync	16	Hot Swap Recognition	24	Hours TMDS-
C1	Red signal			C3	blue signal
C2	Green signal	C5	GND	C4	H-sync

List of cables with connectors

Serial cable 1.8 m (RS232)
USB cable 1.8 m (type A - B)
International power supply 100 ... 240VAC / 5VDC / 3 A
International power supply 100 ... 240VAC / 5VDC / 4 A
VGA cable 1.8m (VGA to DVI-I)
DVI-D cable 1.8 m (DVI-D)
DVI-D splitter cable
HDMI cable 1.8 m
DisplayPort cable 1.8m
Stereo cable with 1.6 m plug (stereo 3.5 mm)
RCA cable 2.5 m (Cinch plug)
TOSLINK cable 1.8 m (F05 plug)
Mini-XLR cable 1.8 m (3-pole)
IEC connection cable for power supply, lockable

Power supply for extension cords

5 VDC	Contact	Signal
	Inside	VCC(+5VDC)
	Outside	GND

KVM Extender :

- Single head devices: max. 800 mA
- Embedded DVI-I devices: max. 1200 mA
- VGA: max. 1700 mA
- Redundancy devices: max. 1050 mA

NET SWITCH 8 8-port network switch (optional) (ProSAFE 8-port Gigabit PoE Smart Managed Switch with 2 Gigabit Fiber Optic SFP article in the manufacturer's QMS, the article is indicated in the operational documentation)

Network protocol and standards

Compatibility

- IEEE 802.3i 10BASE-T
- IEEE 802.3u 100BASE-TX
- IEEE 802.3ab 1000BASE-T
- IEEE 802.3z 1000BASE-X
- IEEE 802.3x full duplex flow control
- IEEE 802.3af (DTE power over MDI)

Power supply

- Power Consumption: 59.3 W maximum

Network ports

- Eight (8) 10/100/1000 Mbps auto-sensing + 2 gigabit SFP uplinks

physical characteristics

- Dimensions: (W x D x H): 236 x 101.6 x 27 mm (9.3 x 4 x 1.07 inches)
- Weight: 0.7 kg (1.55 lbs)

Specifications

- Forwarding modes: store-and-forward.
- Bandwidth: full duplex 20Gbps
- Network latency: less than 4 µs for 64-byte frames in store-and-forward mode for 1000 Mbps to 1000 Mbps transmission.
- Buffer memory: 512 KB internal memory per unit.
- Address database size: 4 KB of Media Access Control (MAC) addresses per system.
- Mean time between failures (MTBF): 157,004 hours (~ 18 years old)

- Acoustic noise: 0 dB

Electromagnetic radiation

- CE mark, commercial
- FCC Part 15 Class B
- VCCI Class B

• C-Tick

Environmental characteristics

- Operating temperature: 0° to 104°F (0° to 40°C)
- Storage temperature: -4° to 158°F (-20° to 70°C)

- Operating humidity: 90% maximum relative humidity, non-condensing
- Storage humidity: 95% maximum relative humidity, non-condensing.
- Operating altitude: 10,000 feet (3,000 m) maximum
- Storage altitude: 10,000 feet (3,000 m) maximum

Status indicators

- Per RJ45 port: connection / activity, speed,

PoE status

- Per SFP port: channel / speed
- Per device: power, PoE max.

Charger

- External power adapter 48 V, 1.25 A adapted to the country of sale

Electromagnetic Immunity

- EN 55024

Safety

- CE mark, commercial
- cUL IEC 950 / EN 60950
- CB

Switch administration

- Static VLAN IEEE 802.1Q (64 groups, static)
- Secure ports
- IEEE 802.1p Class of Service (CoS)
- Port based QoS
- QoS based on DSCP
- DiffServ
- IEEE 802.3ad Link Aggregation (manually or LACP)
- IEEE 802.1w Rapid Spanning Tree

Protocol (RSTP)

- IEEE 802.1s Multiple Spanning Trees

Protocol (MSTP)

- IEEE 802.1ab link layer discovery

Protocol (LLDP)

- LLDP-MED
- SNMP v1, v2c, v3
- RFC 1213 MIB II
- Ethernet interface MIB RFC 1643
- RFC 1493 Bridge MIB
- RMON groups 1, 2, 3, 9
- Auto voice VLAN
- DHCP filtering
- Automatic denial of service (DoS) protection
- HTTP and HTTPS
- Ping and tracing
- Environmentally friendly: energy saving due to cable length (<10 m)
- Power saving when communication is disconnected
- DHCP client RFC 2131
- DHCP filtering
- IEEE 802.1x (MD5) with guest VLAN
- Support for jumbo frames
- Port based security with blocked MAC addresses
- ACL based on MAC and IP
- Storm control for broadcast, multicast and unknown unicast packets.
- IGMP Snooping v1 / v2

- Upload rate limit based on port
- SNTP
- Support for port mirroring (many to one)
- Internet configuration
- Configuration backup / restore
- Password access control
- Support for TACACS+ and RADIUS
- System log
- Ability to update the firmware
- PoE (power headroom 46W)

- port network switch NET SWITCH 4 (optional)

Standards	IEEE 802.3, IEEE 802.3u, IEEE 802.3x, IEEE 802.1p, IEEE 802.3af
Interface	5 RJ45 ports 10/100/1000Mbps with 4 PoE ports (port 1 - port 4) Auto negotiation / AUTO MDI / MDIX
network data	10BASE-T: Category 3, 4, 5 UTP cable (maximum 100m) EIA / TIA-568 100Ω STP (maximum 100m) 100BASE-TX: Category 5 UTP, 5e cable (maximum 100m) EIA / TIA-568 100 Ohm STP (Maximum 100m) 1000BASE-T: Category 5, 5e, 6 or higher UTP cable (Maximum 100m) EIA / TIA-568 100Ω STP (Maximum 100m)
Maximum power consumption	3.95 W (220V/50Hz. No PD connected) 64.44 W (220 V / 50 Hz with 56 W power parts connected)
Maximum heat dissipation	13.48 BTU/h (no PD connected) 219.87 BTU/h (with 56W power parts connected)
PoE ports (RJ45)	Standard: Compliant with 802.3af PoE Ports: Port1 - Port4 Power supply: 56W
switching capacity	10 Gbps
Mac address table	2K
Number of fans	Fanless
External power supply	External Power Adapter (Output: 48V DC / 1.25A)
Dimensions (W x D x H)	3.9 x 3.9 x 1.0 inches (99.8 x 98 x 25 mm)
Certification	FCC, CE, RoHS
Environmental conditions	Operating temperature: 0° C ~ 40° C (32° F ~ 104° F) Storage Temperature: -40° C ~ 70° C (-40° F ~ 158° F)

	Operating Humidity: 10% ~ 90% RH, non-condensing Storage Humidity: 5% ~ 90% RH, non-condensing
--	---

Microphone IP MIC _

A small, versatile clip-on microphone designed specifically for speaking in noisy environments. Its constant firing pattern (cardioid) provides a high level of feedback suppression. The lightweight clasp on this microphone makes it easy to handle. Mini jack plug for ew G3, ew G4 and XS Wireless series bodypack transmitters.

PECULIARITIES

- High speech intelligibility.
- High degree of rejection of feedback
- Easy fastening to clothes.
- Mini jack for bodypack transmitters

Converter principle	pre-polarized condenser microphone
Selection scheme	cardioid
Sensitivity (free field, no load) (1 kHz)	40 mV / Pa
Maximum. SPL at 1 kHz	120 dB
frequency response	50 to 18,000 Hz
Equivalent Noise Level A (DIN IEC 651)	31 dB
Connector	3.5mm jack
Cable length	1.60 m
Supply voltage	7.5V

USB joystick CAM JOYSTICK (AXIS T 8310 article in the manufacturer's QMS, the article is indicated in the operational documentation)

Hall effect joystick with three axes:

X / Y: for pan and tilt,

Z: knob to increase Deviation by square separator, (± 150) for XY,

(± 250) for Z 6 hotkeys with customizable function

USB Type B connectors

dimensions

Height: 109 mm (4.3 inches)

Width: 182 mm (7.2 inches)

Depth: 181 mm (7.1 inches)

Weight 730 g (1.6 lb)

1 USB cable 200 cm (78.7 inches) for PC connection

control **module** (on request)

Analog input channels 8 and 4

Analog output channels 4

Configurable DIO channels (by jumper) 8 and 4

Digital input channels 16,8,6,4

Digital output channels 16,4

Isolation 3k VDC or 2k Vrms

Relay channels 6

RTD channels 6

Thermocouple channels 8

Buttons: Reset button

Digital inputs

Connector Euroblock screw terminal

Counter frequency 250 Hz

Digital filtering time interval : automatic

Dry contact:

On: short circuit to GND

Off: open

I/O mode DI or event counter

Points for COM : 12, 8,6,4 channels

Sensor type: Dry contact

Wet contact (NPN or PNP)

Wet Contact (DI to COM) On: 10 to 30 VDC

Off: 0 to 3 VDC

Digital outputs

Connector Euroblock screw terminal

Current value: 200 mA for all channels per channel and 500 mA for 12 channels per channel

I/O mode DO or pulse output

Overcurrent protection for all channels 2.6A per channel at 25°C and 1.5A for 12 channels per channel at 25°C

Thermal shutdown 175°C (nominal), 150°C (min.)

Surge protection 35 VDC

Pulse output frequency 500Hz (Max.)

Breakdown voltage 500 VAC

Connector Euroblock screw terminal

Contact current rating Resistive load: 5 A at 30 VDC, 250 VAC, 110 VAC

Contact resistance 100 mΩ (max.)

Electrical durability 100,000 operations with 5 A resistive load

Initial insulation resistance 1000 MΩ (min.) at 500 VDC

Mechanical endurance 5,000,000 operations

Pulse output frequency 0.3 Hz at rated load (max.)

Type Form A (N.O.) power relay

Note The ambient humidity must not condense and is between 5 and 95%. The relay may fail when operating in highly condensing conditions below 0°C.

Analog inputs

Accuracy:

±0.1% of full scale at 25°C

±0.3% of full scale at -10 to 60°C

±0.5% of full scale at -40 to 75°C

Built-in resistor for 120 ohm current input

Connector Euroblock screw terminal

Voltage/current in I/O mode

I/O type Differential

Input impedance 10 MΩ (min.)

Input range 0 to 10 VDC

0 to 20 mA

4 to 20 mA

4 to 20 mA (with burnout detection)
 Resolution 16 bit
 Sampling Rate All channels: 12 samples/sec.
 Per channel: 1.5 samples/sec.
 10mA voltage output short circuit protection
 Internal resistor 400 ohm

Note: A 24V external power supply is required if the load exceeds 1000 ohms.

Output range 0 to 10 VDC 4 to 20 mA
 Resolution 12 bits
 Accuracy
 $\pm 0.1\%$ of full scale at 25°C
 $\pm 0.3\%$ of full scale at -10 to 60°C
 $\pm 0.3\%$ of full scale at -40 to 75°C
 Connector Euroblock screw terminal
 Input connection 2- or 3-wire
 Input impedance 625 k Ω (min.)
 Sensor type PT1000 (-200 to 350°C)
 PT50, PT100, PT200, PT500 (-200 to 850°C)
 Resistance type 310, 620, 1250 and 2200 ohm
 Resolution 0.1°C or 0.1 ohm
 Sampling Rate All channels: 12 samples/sec.
 Per channel: 2 samples/sec

Thermocouples

Millivolt Accuracy:

$\pm 0.1\%$ of full scale at 25°C
 $\pm 0.3\%$ of full scale at -10 to 60°C
 $\pm 0.3\%$ of full scale at -40 to 75°C
 TC accuracy types J, T, E, S, B: $\pm 5^\circ\text{C}$
 Types K, R, N: $\pm 8^\circ\text{C}$
 CJC accuracy $\pm 0.5^\circ\text{C}$ at 25°C
 $\pm 1.5^\circ\text{C}$ -40 to 75°C

Input impedance 10 M Ω (min.)

Millivolt type $\pm 19.532\text{ mV}$

$\pm 39.062\text{ mV}$

$\pm 78.126\text{ mV}$

Fault and surge protection: -35 to +35 VDC (power off); -25 to +30 VDC (power on)

Resolution 16 bit

Sampling Rate All channels: 12 samples/sec.

Per channel: 1.5 samples/sec.

Sensor type J, K, T, E, R, S, B, N

Ethernet interface

Ports 10 / 100BaseT (X) (RJ45 connector) 2, 1 MAC address (Ethernet bypass)

Magnetic isolation 1.5 kV (embedded)

Configuration Options: Web Console (HTTP), Windows Utility (ioSearch)

Industrial protocols EtherNet/IP adapter (slave), Modbus TCP server (slave), Moxa AOPC (active tag), MXIO library

BOOTP management, RESTful API, DHCP client, HTTP, IPv4, TCP/IP, UDP, SNMPv1 trap, SNMPv1/v2c

MIB MIB device settings

Security access control list

LED indicators Power, Ready, Port 1, Port 2

Mod TCP bus

Supported functions 1, 2, 3, 4, 5, 6, 15, 16, 23

Mode Server

Maximum. Number of client connections 10

Ethernet/IP

Adapter Mode

Maximum number of scanner connections 9 (read-only), 1 (read/write)

Power parameters

Power connector Euroblock screw terminal

Number of power inputs 1

Input voltage 12 to 36 VDC

Power consumption

16 channels: 110 mA at 24 VDC

8 channels: 155mA at 24VDC

12 channels: 130 mA at 24 VDC

6 channels: 188 mA at 24 VDC

4 channels: 139 mA at 24 VDC

physical characteristics

Housing Plastic

Dimensions 27.8 x 124 x 84 mm (1.09 x 4.88 x 3.31 inches)

Weight 200 g (0.44 lb)

Installation DIN rail mounting, Wall mounting

I/O cable connection, 16 to 26 AWG

Power cord, 12 to 24 AWG

Environmental conditions

operating temperatures: -10 to 60°C (14 to 140°F)

Storage Temperature -40 to 85°C (-40 to 185°F)

CAMERA IP COMPACT video camera (optional)

Use with the following browsers:

- Chrome™ (recommended), Firefox® , Edge® or Opera® with Windows®
- Chrome™ (recommended) or Safari ® with OS X® _
- Chrome™ or Firefox ® with other operating systems.

Microphone 3.5mm

Audio extension cable A , 5 m (16 ft)

16 image frames per second. Full video at 30 (NTSC) or 25 (PAL) frames per second.

H.264 or MPEG -4 Part 10/ AVC

Output 12V DC


I/O connector

Function	Contact	Notes	Characteristics
0 V DC (-)	GND	DC ground	0 VDC
DC output	12 V	Can be used to power auxiliary equipment. Note. This pin can only be used for power output.	12 VDC Max load = 50mA
Customizable (input or output)	I/O 1-4	Digital Input - Connect to the GND pin to enable or leave unconnected to disable.	0 to max. 30 V DC

		Digital output - internally connected to pin 1 (DC ground) when active and not connected when inactive. When used with an inductive load, such as a relay, a diode should be connected in parallel with the load to protect against power surges.	0 to 30 VDC max, open drain, 100 mA
--	--	--	-------------------------------------

Power connector

2-pin terminal block for DC input. Use a Limited Power Supply (LPS) compliant with Safety Extra Low Voltage (SELV) either with rated output power limited to $\leq 100W$ or rated output current limited to $\leq 5A$.

CAMERA IP HD OUTDOOR video camera (optional)

Identification marks acc. to Directive 2014/34/EU:

II 2G (zone 1 and 2)

II 2D (zone 21 and 22)

Explosion protection (gas): Ex db IIC T6 Gb

Explosion protection (dust): Ex tb IIIC T80 °C Db

Protection class: IP 68 (IEC / EN 60529)

Transport / storage temperature: 0 °C... + 40 °C

Ambient temperature (EX) 1: -50°C...+60°C

Power input 24 VDC:

Permissible temperature range: -0 °C < Tamb. < +60°C

Power supply: PoE, IEEE 802.3at class 4

Reference voltage: 48 V DC (44 ... 54 V DC)

Maximum power consumption: 17.1W

Typical power consumption: 11W

Connection cable Ex-d - Ex-e (ASKDP03-T)

Description: Communication and power supply of the camera module (according to DIN EN 60079-14)

Sheath color: green (GN), similar to RAL3001

Outside Diameter: 15.50±0.6mm

Bend radius: 15 x Yes when installed and 10 x Yes after moving

Data line: 4 x 2 x AWG23 / 1 CAT.6a

Performance Elements: 3G1.5 (BK-BU-GN/YE)

	Camera (Ex-d)	Terminal block (Ex-e)
Permissible ambient temperature:	0°C...+60°C (for PoE+ power supply) -50°C...+60°C (when powered by 24V DC)	-60°C...+55°C
Protection class according to EN 60529 / IEC 529:	IP68 (Test conditions: 24h / 3m, w.st 5 °C °)	IP66
Body material:	stainless steel, mat. number 1.4404 LEXAN	polyester resin
Weight:	20 kg	about 1 kg

Dimensions:	D195mm x 378mm	145mm x 145mm x 71mm
-------------	----------------	----------------------

CAMERA IP HD PRO video camera (optional)

Use with the following browsers:

- Chrome™ (recommended), Firefox®, Edge® or Opera® with Windows®
- Chrome™ (recommended) or Safari ® with OS X® _
- Chrome™ or Firefox ® with other operating systems.

Microphone 3.5mm

Audio extension cable A , 5 m (16 ft)

16 image frames per second. Full video at 30 (NTSC) or 25 (PAL) frames per second.

H.264 or MPEG -4 Part 10/ AVC


Output 12VDC

I/O connector

Function	Contact	Notes	Characteristics
0 V DC (-)	GND	DC ground	0 V DC
DC output	12 V	Can be used to power auxiliary equipment. Note. This pin can only be used for power output.	12 V DC Max load = 50 mA
Customizable (input or output)	I/O 1-4	Digital Input - Connect to the GND pin to enable or leave unconnected to disable.	0 to max. 30 VDC
		Digital output - internally connected to pin 1 (DC ground) when active and not connected when inactive. When used with an inductive load, such as a relay, a diode should be connected in parallel with the load to protect against power surges.	0 to 30 VDC max, open drain, 100 mA

Power connector

2-pin terminal block for DC input. Use a Limited Power Supply (LPS) compliant with Safety Extra Low Voltage (SELV) either with rated output power limited to $\leq 100W$ or rated output current limited to $\leq 5A$.

Video camera CAMERA IP HD PLUS

CAMERA IP HD PLUS	Power supply	Dedicated power adapter	CAMERA IP PLUS POWER
	Audio	External microphone connected directly to the camera	IP MIC

	Control Protocol	VAPIX (axle)	-
--	------------------	--------------	---

BNC connector

The BNC connector is used to connect broadcast equipment using SDI. Connect a 75 ohm SDI coaxial cable.

The HDMI connector is used to connect to a display or monitor.

I/O connector

Use with external devices in conjunction with, for example, tamper alarm, motion detection, event trigger, frame recording and alarm notification. In addition to the 0 VDC reference point and power (DC output), the I/O connector provides an interface for:

Digital output - for connecting external devices such as relays and LEDs. Connected devices can be activated using VAPIX® Application Programming Interface, output buttons on live view page or via action rule. The output will show as active (displayed under System Settings > Ports and Devices) if the alarm device is activated.

- Digital input - an alarm input for connecting devices that can switch between open and closed circuits, for example: PIR, door/window contacts, glass break sensors, etc. When a signal is received, the state changes and the input becomes active (displayed under System Settings > Ports and Devices)

Power connector

2-pin terminal block for power input. Use a Limited Power Supply (LPS) compliant with Safety Extra Low Voltage (SELV) either with rated output power limited to $\leq 100W$ or rated output current limited to $\leq 5A$.

Audio jack

The product has the following audio connectors:

- Audio input - 3.5 mm stereo microphone or stereo line level input.
- Audio output - 3.5 mm audio output (line level) that can be connected to a public address system (PA) or powered speaker with built-in amplifier. You can also attach a pair of headphones. You must use a stereo jack for audio output.
- Left - 3-pin XLR connector for balanced audio input. Use the left jack for mono.
- Right - 3-pin XLR connector for balanced audio input.

Has a slot for an SD card.

Name	Characteristic	Data
Camera	Image sensor	r 1/2.8" progressive scan RGB CMOS
	Lens	Autofocus. P-iris control, 4.3 - 129mm, F1.6-4.7 Horizontal viewing angle: 63.7° - 2.3° Vertical viewing angle: 38.5° - 1.3°
	Minimum lighting	Color: 1.0 Lux at 30 IRE, F1.6 B/W: 0.03 lux at 30 IRE, F1.6
	holding time	1/10000 s to 1 s
	Pan / Tilt / Zoom	Rotation: $\pm 170^\circ$, 0.2-100°/s Tilt: -20° - $+90^\circ$, 0.2 - 90°/s Accuracy: $\pm 0.2^\circ$ Resolution: 0.02° Zoom: 30x optical. 12x digital, 360x total

	Pan/Tilt/Zoom functions	256 preset positions, control queue, on-screen direction indicator, adjustable zoom and speed
Video	Video compression	H.264 High, Baseline and Main profile (MPEG-4 Part 10 / AVC) Motion JPEG
	Resolutions	1920x1080 (HDTV 1080p) to 320x180 HDMI / SDI: up to HDTV 1080p
	Personnel	Up to 50/60 fps (50/60 Hz)
	Video broadcast	Multiple streams in H.264 and Motion JPEG. H.264 and Motion JPEG: 1 stream at full resolution and frame rate. More streams if they are identical or limited in frame rate or resolution. Adjustable frame rate and bandwidth VBR/MBR H.264
	Image settings	120dB WDR - dynamic capture, compression, color, brightness, sharpness, white balance, exposure control, exposure zones, backlight compensation, fogging, noise reduction, electronic image stabilization, text and image overlay (via Ethernet) Rotation: 0°, 180°
Audio	Streaming Audio	Two-way, stereo
	Audio compression	SDI: AES3 24bit, 48kHz HDMI: LPCM Net: AAC LC 8 to 48 kHz, 8 to 128 kbps G.711 PCM 8 kHz 64 kbps G.726 ADPCM 8 kHz 24 to 32 kbps Opus 8/16/48 kHz 8 to 128 kbps Configurable baud rate
	Audio input	Balanced mono XLR input (left/right) Balanced external microphone (mono) 48V phantom microphone power Balanced line level (mono) Line input impedance: >7.1 kΩ Maximum input level: 4.4 Vrms Bandwidth: 20Hz - 20kHz, can be sample rate limited


		THD+N: <0.03% Signal-to-noise ratio: >85dB at 0dB gain, >78dB at 30dB gain Unbalanced 3.5mm stereo input Unbalanced external microphone (stereo) Microphone power supply 5 V via 2.2 kOhm Unbalanced line (stereo) Line input impedance: >1.7 kΩ Maximum input level: 2.2 Vrms THD+N: <0.03% Signal-to-noise ratio: >87dB at 0dB gain, >83dB at 30dB gain Bandwidth: 20Hz - 20kHz, can be sample rate limited
	Audio output	Unbalanced 3.5mm stereo output Output impedance: <100 ohm, short circuit protection Signal to noise ratio: >87dB Bandwidth: 20Hz - 20kHz, can be sample rate limited THD+N: <0.03% at 10kΩ load Maximum output level: >0.707 Vrms SDI THD+N: <0.03% Bandwidth: 20Hz - 20kHz Signal-to-noise ratio: >87dB at 0dB gain, >83dB at 30dB gain HDMI THD+N: <0.03% Bandwidth: 20Hz - 20kHz Signal-to-noise ratio: >87dB at 0dB gain, >83dB at 30dB gain
Net	Safety	Password protection, IP address filtering, HTTPS encryption, IEEE 802.1X network access
	Supported protocols	IPv4/v6, HTTP, HTTPS*, SSL/TLS*, QoS Layer 3 DiffServ, FTP, SMTP, Bonjour, UPnP™, SNMP v1/v2c/v3 (MIB-II), DNS, DynDNS, NTP, RTSP, RTP, TCP, UDP, IGMP, RTCP, ICMP, DHCP, ARP, SOCKS, CIFS/SMB, SSH, HDMI 1.4, SDI
System integration	Built-in installers	Pixel counter
	Streaming data	Event Data

	Event triggers	Analytics, Edge storage events, external input
	Event actions	File upload via FTP, SFTP, HTTP, HTTPS, network share and email notification. Email notification, HTTP, HTTPS and TCP. Sending a video clip, sending an SNMP trap. Record video and audio to peripheral storage. Video buffering before and after an alarm. External output activation Day/Night Mode Audio Clip Playback
	Video access	Live view from camera Write video to file (ASF) Custom HTML pages Windows 7, Windows Vista, Windows XP, Windows Server 2008, Windows Server 2012 and Windows 8 DirectX 9c or higher
Are common	Frame	Plastic cover ASA black and white. Colour: white NCS S 1002-B.
	Memory	512 MB Flash, 1 GB RAM Real time clock with battery backup
	Power	8-28VDC (including 12V power supply) typ. 12W, max. 21 W
	Connectors	RJ45 10BASE - T / 100BASE - TX Terminal block for 2 configurable alarm inputs/outputs 3.5mm mic/line input jack 3.5mm stereo line out jack XLR-3 (left + right) mic/line in (with 48V phantom power) HDMI Type A BNC to SDI, SMPTE 292-1 (1.5Gbps) or SMPTE 424 (3Gbps) depending on capture Mode Terminal block for DC power supply
	Edge Storage	SD/SDHC/SDXC slot supporting up to 64GB memory card (card not included) Support for writing to dedicated network storage (NAS)

	terms of Use	Temperature: 0°C to 40°C (32°F to 104°F) Humidity 10-85% (non-condensing)
	Storage conditions	Temperature: -20 to 60°C (-4 to 140°F) Humidity: 10-85% (non-condensing)
	Statements	EN 55022 class A, EN 55024, EN 61000-6-1, EN 61000-6-2, EN 61000-3-2, EN 61000-3-3 FCC Part 15 Subpart B Class A ICES-003 Class A VCCI Class A KSS KN 22 Class A, KN 24 C-tick AS/NZS CISPR 22 Class A IEC/EN/UL 60950-1 EN 50581
	dimensions	180 x 136 x 136 mm (7.1 x 5.4 x 5.4 inches)
	Weight	1.49 kg (3.3 lb)
	Languages	German, French, Spanish, Italian, Russian, Simplified Chinese, Japanese, Korean, Portuguese
	Guarantee	3 years

Components of the Micromed system:

Component	Characteristics
Bracket ARM GM LTM	Two-knee, knees - 71 and 92 cm, 4.7 kg
Horizontal bracket ARM LCD HORIZ	61 x 63 cm, 4.2 kg.
Vertical bracket ARM LCD SW	Height - 60 cm, weight - 2 kg.
Adjustable BAG BELT HP	Length - 120 cm, weight - 300 gr.
Bag LTM HP	15 x 5 x 10 cm, weight - 110 gr.
Rack with base BASE STAND (optional)	Assembled - 1 m, extended - 1.6 m, whole - 10 kg
CAM STAND	In the assembled state - 70 cm, in the extended state - 1 m, weight - 3.5 kg
Wheeled base LED STAND PHOTIC	5 wheels, wheel diameter 10 cm
Bracket LED SUPP PHOTIC	30 x 20 x 10 cm, weight 900 gr.

Stand STAND POLE NEW on wheels with weighting (optional)	5 wheels, wheel diameter 10 cm, weight 5 kg
ADAPT DVI HDMI adapter (optional)	DVI to HDMI adapter
Tripod stand VSTAND	Height - 1.2 m, extended - 1.8 m, weight - 5 kg.
Wall bracket VSUPP	40 cm length, weight - 2 kg.
VSUPP IP metal support for mounting IP video cameras	Diameter - 30 cm, weight - 340 gr.
ADAPT BNC RCA adapter (optional)	BNC to RCA adapter
Adapter CB IP HD PLUS EXTENDER	As in the case of HBC BT/RJ45- IIBC HT/RJ45 consists of adapters, but the manufacturer has one article CB IP HD PLUS EXTENDER for this set of two. The function is extension. Block 1 - Com Port - RJ 45 block 2 - RJ 45- DIN 8 Pin . Photo of CB IP HD PLUS EXTENDER attached
Trolley CAR STAR	Height - 70 cm, Width - 50 cm, Weight - 12 kg.
Trolley CAR PLUS	Height 120 cm, Width - 50 cm, weight - 15 kg
Trolley CAR TORO	Height - 60 cm, Width 50 cm, weight - 9kg
Trolley CAR TORO SMALL	Height - 45 cm, Width 40 cm, weight - 7kg
NET ISOLAT Network Isolator for 10/100 Mbps Networks	The manual is attached, in the characteristics it should be noted that it only works with 10/100 Mbps
Professional network switch NET SWITCH PRO (optional)	Manual attached
Forceps PLIERS TP	The tweezers are plastic.  Length 150 mm +/-10%
Patient button with 5m cable and RCA connector (optional)	Plastic button with 5 m copper cable and RCA connector
Illuminated patient button with 2 m cable and LEMO connector (optional)	Plastic, with backlight and LEMO connector (quick-detachable connector with 4 pins) for connection to the Wall socket plug . Cable length - 1.8 mm.

Power supply for BUTTON BOX (optional)	Input 100-240 V , 50-60Hz. 0.3A, output - +9V DC, 1.1 A
wall panel EXPRESS WALL CONNECTOR NT /DE with	Width-25cm. height 20cm. weight- 3kg
HBC BT/RJ45 extension adapter	HBC BT/RJ45 is used in conjunction with HBC HT/RJ45, extension function. If you need an extension cord, then this adapter is connected with one end (4 wires) to the terminals of the Terminal Box, and with the other end to the patch cord. The second end of the patch cord is connected to the HBC HT/RJ45, the other end of the HBC HT/RJ45 is connected to the computer.
HBC HT/RJ45 extension adapter	
EXTERNAL BATTERY (optional)	Battery type - LI - ION , charging time - 3 hours, operating time on a single charge - no more than 24 hours, capacity - 3600 m / Ah, voltage - 5 V, overall dimensions - 10 length, diameter - 3 cm, weight - 300 gr, the battery is removable.
EXTERNAL CHARGER	Consumed power - 50 watts, weight - 300 gr, 10 x 7 x 7 cm, connector type - USB - micro USB
Bracket ARM GM MORPH	Two-knee, knees - 68 and 87 cm, 4.2 kg
BAG MORPH bag	20 x 10 x 7 cm, 200 gr.

BQ NET POE 2 interface module

Output Trigger Specifications (RED CONNECTOR)

One trigger output on RCA jack

Amplitude: (0-5V) positive pulses

duration: 30 μ s

The output is isolated by a pulse transformer.

It is possible to change the pulse polarity by inverting the cable pole.

Trigger Input Specifications (BLACK CONNECTOR)

One trigger input on RCA jack

Amplitude: (0-5V)

Rising edge trigger if trigger duration is less than 25 μ s.

Maximum trigger frequency in fs/2 (1 trigger for every 2 max samples).

There is no minimum trigger duration because the edge is recognized (recommended 10 μ s) Isolated trigger with pulse transformer. It is possible to accept a falling edge impulse by inverting the cable pole.

BQ USB interface module

Output Trigger Specifications (RED CONNECTOR)

One trigger output on RCA jack

One USB connector

Amplitude: (0-5V) positive pulses

duration: 30 μ s

The output is isolated by a pulse transformer.

It is possible to change the pulse polarity by inverting the cable pole.

Trigger Input Specifications (BLACK CONNECTOR)

One trigger input on RCA jack

Amplitude: (0-5V)

Rising edge trigger if trigger duration is less than 25 μ s.

Maximum trigger frequency in fs/2 (1 trigger for every 2 max samples).

There is no minimum trigger duration because the edge is recognized (recommended 10 μ s) Isolated trigger with pulse transformer. It is possible to accept a falling edge impulse by inverting the cable pole.

BQ interface module NET POE PLUS

Output Trigger Specifications (RED CONNECTOR)

One trigger output on RCA jack

Amplitude: (0-5V) positive pulses

duration: 30 μ s

The output is isolated by a pulse transformer.

It is possible to change the pulse polarity by inverting the cable pole.

Trigger Input Specifications (BLACK CONNECTOR)

One trigger input on RCA jack

Amplitude: (0-5V)

Rising edge trigger if trigger duration is less than 25 μ s.

Maximum trigger frequency in fs/2 (1 trigger for every 2 max samples).

There is no minimum trigger duration because the edge is recognized (recommended 10 μ s) Isolated trigger with pulse transformer. It is possible to accept a falling edge impulse by inverting the cable pole.

BQ NET POE M interface module

Output Trigger Specifications (RED CONNECTOR)

One trigger output on RCA jack

Amplitude: (0-5V) positive pulses

duration: 30 μ s

The output is isolated by a pulse transformer.

It is possible to change the pulse polarity by inverting the cable pole.

Trigger Input Specifications (BLACK CONNECTOR)

One trigger input on RCA jack

Amplitude: (0-5V)

Rising edge trigger if trigger duration is less than 25 μ s.

Maximum trigger frequency in fs/2 (1 trigger for every 2 max samples).

There is no minimum trigger duration because the edge is recognized (recommended 10 μ s) Isolated trigger with pulse transformer. It is possible to accept a falling edge impulse by inverting the cable pole.

BQ power supply NET POWER POE

Power supply: 37-57VDC from POE

Power Consumption: <6.49-12.95W

Maximum output power: 5W

Trigger isolation: 3000V with transformer

Amplifier data isolation: 2500V with optocouplers

HDMI EXTENDER device

Resolution Supported: 480i@60Hz, 480p@60Hz, 576i@50Hz, 576p@50Hz, 720p@50/60Hz, 1080i@50/60Hz, 1080p@50/60Hz

- Mains cable: Cat.6 / 6A / 7 excellent quality copper
- Supported transmission distance: 60m
- Maximum supported video resolution: 1080P
- Supported Bandwidth: 10.2Gbps

- Supports 8-bit deep color
- Supports PCM stereo audio
- HDCP 1.3 compliant
- Consumption: TX: 5W; RX: 5W
- Power supply: DC5V1A * 2 pcs.
- Product dimensions: LxWxH: 66x56x22 mm.

ANYWHERE 2 device

dimensions

Length: 2.38 inches (6.04 cm)

Width: 3.9 inches (10cm)

Height: 1.0 inches (2.54 cm)

Weight: 5 oz. (142 g)

Environmental conditions

Operating Temperature: 32°F to 131°F (0°C to 55°C)

Relative humidity: 0% to 95% (non-condensing)

Anywhere 2 uses a 120/230VAC, 50/60Hz power adapter that supplies 5VDC to the device.

We recommend that you use only the supplied AC adapter with Anywhere 2.

Network Interface Features

Standards: IEEE 802.3, 802.3i (10Base-T), 802.3u (100Base-TX), 802.3x (full duplex and flow control), HP Auto-MDIX (automatic detection of straight or crossover cable)

Physical layer: 10/100 Mbps half duplex or full duplex, auto-negotiation of speed and duplex .

Ethernet connector: RJ-45

The device has 2 USB ports (standard type A sockets). Downstream ports support low, full, and high speed downstream devices.

RAM: 64 MB

ANYWHERE 5 device

dimensions

Length: 4.35 inches (11.05 cm)

Width: 7.20 cm (18.29 inches)

Height: 1.03 inches (2.61 cm)

Weight: 10.00 oz. (283.5 g)

Environmental conditions

Operating Temperature: 32°F to 131°F (0°C to 55°C)

Relative humidity: 0% to 95% (non-condensing)

Anywhere 5 uses a 120/230VAC, 50/60Hz power adapter that supplies 5VDC to the device.

We recommend that you use only the supplied power adapter with Anywhere 5.

The device has 5 USB ports (standard type A sockets). Downstream ports support low, full, and high speed downstream devices.

Memory: 128MB RAM

Standards: IEEE 802.3, 802.3i (10Base-T), 802.3u (100Base-TX), 802.3x (full duplex and flow control), HP Auto-MDIX (automatic detection of straight or crossover cable)

Physical layer: 10/100 Mbps half duplex or full duplex, with auto-negotiation of speed and duplex.

Ethernet connector: RJ-45

NET ISOLAT Network Isolator for 10/100 Mbps Networks

- High performance Gigabit Ethernet
- Dielectric strength 5.0 kV AC
- Dielectric strength 8.5 kV DC
- Complies with IEC 60601-1

- UL recognized component
- ISO 11801 Class D Ethernet performance achievable
- Extremely low insertion loss; a total cable length of 100 meters is thus achievable.
- ESD protection: transient voltage suppression on all signal lines
- Suitable for devices with supply voltage up to 400 VAC.
- RoHS Compliant

Annex E

(reference)

Table of software options. Description of the software option and correlation of options with screenshots (images of the working view of the software option in working conditions)

No.	Software option	vendor code	Screenshot #	Option description
1	Option for software recording and viewing of EEG, 2-channel CFM version	SW ACQ CFM2	46	Allows registration, viewing and analysis of two channels (C3 and C4) of EEG recording when monitoring the cerebral function of the brain
2	Option for software recording and viewing of EEG, 8-channel CFM version	SW ACQ CFM8	45	Allows registration, viewing and analysis of eight channels (Fp 1, Fp 2, T 3, C3, C4, T 4, O 1, O 2) of EEG recordings when monitoring the cerebral function of the brain
3	Option for software recording and viewing of EEG, version: CLINIC	SW ACQ CLINIC	1.43	Allows registration and recording of EEG on 32-channel amplifiers. This full version allows you to register, record and analyze both low frequency (sampling frequency -64, 128, 256 Hz) and HFO EEG (high frequency, sampling frequency 512, 1024, 2048, 4096, 8192, 16384, 32768 Hz) and contains an extended set of tools (filters, markers, selections, etc.) for working with data
4	Option for software recording and viewing of EEG, version: STUDIO	SW ACQ STUDIO	1.43	Allows registration and recording of EEG on 32-channel amplifiers. This is a reduced version that allows you to register, record and analyze low-frequency EEG (sampling frequency -64, 128, 256 Hz), and contains a basic set of tools for working with data.
5	Option for software registration and viewing of EEG, version: MULTICH	SW ACQ MULTICH	1.43	Allows you to register and record EEG using a combination of one, two or four 64-channel amplifiers. This full version allows you to register, record and analyze both low frequency (sampling

				frequency -64, 128, 256 Hz) and HFO EEG (high frequency, sampling frequency 512, 1024, 2048, 4096, 8192, 16384, 32768 Hz) and contains an extended set of tools (filters, markers, selections, etc.) for working with data
6	Software option for recording and viewing digital video-EEG up to 64 channels, CLINIC version	SW ACQ DV64	1.17,44	Allows registration, recording and viewing of digital video in HD (1280x720) format during EEG registration up to 64 channels. The full version, which contains the ability to view videos with increased and decreased playback speeds, and also has a digital zoom.
7	Software option for recording and viewing digital video-EEG, CLINIC version	SW ACQ DVCL	1.17	Allows registration, recording and viewing of digital video in HD format (1280x720) during EEG registration up to 32 channels. The full version, which contains the ability to view videos with increased and decreased playback speeds, and also has a digital zoom.
8	Software option for recording and viewing digital video-EEG in Full HD format	SW ACQ DVHD	1.17	Allows registration, recording and viewing of digital video in FullHD format (1920 x 1080) when registering EEG up to 256 channels. The full version, which contains the ability to view videos with increased and decreased playback speeds, and also has a digital zoom. Does not have the ability to write video from cameras connected via TCP / IP protocol .
9	Software option for recording and viewing digital IP video-EEG over a network, CLINIC version	SW ACQ DVIP	1.42	Allows registration, recording and viewing of digital video in FullHD (1920 x 1080) format from cameras connected via TCP / IP protocol during EEG registration up to 32 channels. The full version, which contains the ability to view videos with increased and decreased playback speeds, and also has a digital zoom.
10	Software option for recording and viewing digital IP video-EEG over a network of up to 64 channels, CLINIC version	SW ACQ DVIP64	1.42	Allows registration, recording and viewing of digital video in HD format (1280 x 720) from cameras connected via IP protocol during EEG registration up to 64 channels. The full version, which contains the ability to view videos

				with increased and decreased playback speeds, and also has a digital zoom.
eleven	Software option for recording and viewing EEG and video over the network in real time	SW ACQ DVNET	4,5,6	Allows access from the viewing station to the EEG and video registration and recording station via a local network with the ability to remotely control the process of recording and recording both EEG signals and digital video.
12	Software option for EEG analysis in "on-line" and "off-line" modes, version: full	SW ANALYSER ACQ	2,3,35,36	A specialized tool for analyzing EEG and polygraphic signals both in registration and recording mode and in viewing mode. It provides quantitative data analysis, both in the frequency and time domain, showing the evolution of any of the calculated parameters over time. The full version includes DSA (spectral density table) and CFM (cerebral brain function monitoring) versions.
13	Software option for EEG analysis in "on-line" and "off-line" modes, version: CFM	SW ANALYSER ACQ CFM	2	A specialized tool for analyzing EEG and polygraphic signals both in registration and recording mode and in viewing mode. It provides quantitative data analysis, both in the frequency and time domain, showing the evolution of any of the calculated parameters over time. CFM version (cerebral brain function monitoring).
14	Software option for EEG analysis in "on-line" and "off-line" modes, version: DSA.	SW ANALYSER ACQ DSA	3	A specialized tool for analyzing EEG and polygraphic signals both in registration and recording mode and in viewing mode. It provides quantitative data analysis, both in the frequency and time domain, showing the evolution of any of the calculated parameters over time. DSA version (Spectral Density Table).
15	Software option for transmitting EEG data via TCP/IP	SW EEG BCI		On-line procedure for sending EEG data via TCP/IP to a remote server with BCI interface (brain-computer interface)
16	Software option "ERP-manager"	SW EEG ERP	50	Ability to display a digital trigger code with or without a vertical bar

				the ability to control a mixed trigger input, including one digital and one analog trigger
17	Software option "Management of the outpatient registrar"	SW HOLTER	37	Allows direct viewing of registration in the EEG registration mode, the ability to read data from a card reader, the ability to load data from a memory card without removing it.
18	EEG software viewing option, version: CLINIC	SW REV CLINIC	38.39	Allows viewing, editing and analysis of recorded EEG signals up to 32 channels. This full version contains an extended set of tools (filters, markers, selections, etc.) for working with data.
19	EEG software viewing option, version: STUDIO	SW REV STUDIO	38.39	Allows viewing, editing and analysis of recorded EEG signals up to 32 channels. This is a reduced version that contains a basic set of tools for working with data.
20	EEG software viewing option, version: CLINIC MULTICH	SW REV MULTICH	38.39	Allows viewing, editing and analysis of recorded EEG signals up to 256 channels. This full version contains an extended set of tools (filters, markers, selections, etc.) for working with data.
21	EEG software view option, version: EASY	SW REV EASY	38.39	Allows viewing, editing and analysis of recorded EEG signals up to 24 channels. This is a reduced version that contains a basic set of tools for working with data.
22	Software option for viewing EEG in "off-line" mode, version: full	SW ANALYSER REV	40.41	A specialized tool for the analysis of EEG signals and polygraphic signals in view mode only. It provides quantitative data analysis, both in the frequency and time domain, showing the evolution of any of the calculated parameters over time. The full version includes DSA (spectral density table) and CFM (cerebral brain function monitoring) versions.
23	Software option for viewing EEG in "off-line" mode, version: CFM	SW ANALYSER REV CFM	40.41	A specialized tool for the analysis of EEG signals and polygraphic signals in view mode only. It provides quantitative data

				analysis, both in the frequency and time domain, showing the evolution of any of the calculated parameters over time. CFM version (cerebral function monitoring)
24	Software option for viewing EEG in "off-line" mode, version: DSA	SW ANALYSER REV DSA	40.41	A specialized tool for the analysis of EEG signals and polygraphic signals in view mode only. It provides quantitative data analysis, both in the frequency and time domain, showing the evolution of any of the calculated parameters over time. DSA version (Spectral Density Table).
25	Software option for connecting to databases of external applications	SW DB LINK	block diagram 32,33,34	SysstemPlusEvolution program to work with a plug-in that is installed on the server and communicates with the patient database of a medical institution.
26	EEG software option "Event identification"	SW EEG EVENTS	27.28	Allows the user to use a waveform detector with controllable search parameters and scoping. The detector has the ability to select the waveform for automatic search.
27	EEG software option "Analysis of high-frequency oscillations"	SW EEG HFO	thirty	Allows the user to view and analyze previously recorded high frequency oscillations (sampling rate above 256 Hz).
28	EEG software option "Brain Mapping", version: STUDIO	SW EEG MAPS	24,51,52	Enables the user to use the brain mapping tool. This is a frequency and amplitude analysis of the obtained EEG data with visualization on a map of the brain. The basic version of the software option containing the initial set of tools, filters and functions.
29	EEG software option "Brain Mapping", version: CLINIC	SW EEG MAPCL	24, 51.52	Enables the user to use the brain mapping tool. This is a frequency and amplitude analysis of the obtained EEG data with visualization on a map of the brain. Full version of the software option containing an extended set of tools, filters and functions.
thirty	EEG software option "Data exchange with MRI"	SW EEG MRI	22,23,25	Allows you to automatically remove artifacts from the received EEG data using various methods. The software option can be used both during registration and

				recording, and in the mode of viewing the received EEG data.
31	EEG software option "HOLBERG SCORE interface"	SW EEG SCORE	31	Allows you to transfer the recorded EEG data to a third-party automatic analysis program. Does not contain software, just creates a link to a third party program to create a link between them.
32	EEG software option "Audio playback"	SW EEG SOUND	20,21,26	Allows playback of recorded EEG data in audio format. It is also possible to reproduce both all the obtained "curves" and individually.
33	Software EEG option "Averaging in "off-line" mode"	SW EEG TRIG	19	Gives the user the ability to edit digital and analog triggers (adding, deleting any kind of triggers and synchronizing them with the received EEG data).
34	Software option for viewing digital video-EEG, version: CLINIC	SW REV DVCL	17.18	Allows you to view digital video in HD (1280x720) format of saved EEG data up to 32 channels. The full version, which contains the ability to view videos with increased and decreased playback speeds, and also has a digital zoom.
35	Software option for viewing digital video-EEG, version: STUDIO	SW REV DVST	17.18	Allows you to view digital video in HD (1280x720) format of saved EEG data up to 32 channels. The full version, which contains the ability to view videos with increased and decreased playback speeds, and also has a digital zoom.
36	Software option for viewing digital video-EEG in FULL HD format	SW REV DVHD	17.18	Allows you to view digital video in FullHD format (1920 x 1080) of stored EEG data up to 256 channels. The full version, which contains the ability to view videos with increased and decreased playback speeds, and also has a digital zoom.
37	Software option for viewing digital video-EEG in multichannel mode, version: CLINIC	SW REV DVMULTICH	17.18	Allows you to view digital video in HD (1280x720) format of stored EEG data up to 256 channels. The full version, which contains the ability to view videos with increased and decreased playback speeds, and also has a digital zoom.

38	Software option for advanced management of the archive server	SW S98 SERVER	14	Allows the user to set up automatic archiving of decrypted EEG data.
39	Program statistics option	SW S98 STAT	15.16	Allows an application to work with a database created in a Microsoft application SQLServer. Thanks to this option, it becomes possible to maintain detailed statistics of patient studies.
40	Software option for the sleep program interface	SW S98 SW	7	Allows you to transfer the recorded EEG data to a third-party program for automatic analysis. Does not contain software, just creates a link to a third party program to create a link between them.
41	Software option for background data transfer	SW S98 TRANSFER	11.12	This option allows the user to configure the automatic transfer of recorded EEG data between recording stations and the server.
42	Software option for exporting EEG data to the Viewer	SW S98 VIEWER	13	Allows you to export saved studies along with a program for viewing exported EEG data.
43	Software option to enable HIS interface	SW XML INTERFACE	block diagram 32,33,34	SysystemPlusEvolution program to work with a plug-in that is installed on the server and communicates with the patient database of a medical institution.
44	Software option for analyzing imported curves in EDF format	RT OPT EDF	10	Allows the application to export study files in EDF format for subsequent import into a third-party automatic sleep analysis program.
45	Software option for analyzing curves recorded on EEG head units for long-term monitoring SD LTM 32 EXPRESS (32 channels)	RT OPT SD LTM EXP32	7.47	Allows you to export saved EEG data recorded on 32-channel LTM series amplifiers EXPRESS without formatting into a third-party program for automatic sleep analysis.
46	Software option for analyzing curves recorded on EEG head units for long-term monitoring SD LTM 64 EXPRESS (64 channels)	RT OPT SD LTM EXP64	7.47	Allows you to export saved EEG data recorded on 64-channel LTM series amplifiers EXPRESS without formatting into a third-party program for automatic sleep analysis.
47	Software option for analyzing curves recorded on EEG Head Units for long-term monitoring 32 channels	RT OPT SD LTM32	7.47	Allows you to export saved EEG data recorded on 32-channel amplifiers without formatting to a third-party program for automatic sleep analysis

48	Software option for analyzing curves recorded on EEG Head Units for long-term monitoring 64 channels	RT OPT SD LTM64	7.47	Allows you to export saved EEG data recorded on 64-channel amplifiers without formatting to a third-party program for automatic sleep analysis.
49	Software option for the analysis of curves recorded on EEG Head Units of the SD PLUS series	RT OPT SD PLUS	7.47	Allows you to export saved EEG data recorded on SD series amplifiers PLUS without formatting into a third-party program for automatic sleep analysis.
50	Software option "Centralized control of video cameras" on a USB key	SW CAMERA CONTROL	9	Allows the user to remotely control video cameras from which recording and recording of video synchronized with EEG data is performed.
51	Software option "Event management" (on demand) on a USB key	SW EVENT MANAGER	8	Thanks to this software option, it is possible to receive and display various events (generated by a patient button, an external button or an EEG analyzer) coming from different data acquisition devices.
52	Software option for recording and viewing digital video-EEG over the network, CLINIC version		1.17	Allows registration, recording and viewing of digital video in HD (1280x720) format during EEG registration from cameras connected via TCP / IP protocol , up to 32 channels. The full version, which contains the ability to view videos with increased and decreased playback speeds, and also has a digital zoom.
53	Software option for recording and viewing digital video-EEG over a network of up to 64 channels, CLINIC version		1.17	Allows registration, recording and viewing of digital video in HD (1280 x 720) format during EEG registration from cameras connected via TCP / IP protocol up to 64- x channels. The full version, which contains the ability to view videos with increased and decreased playback speeds, and also has a digital zoom.
54	Software option for recording and viewing digital video-EEG up to 128 channels, version: CLINIC	SW ACQ DV128	1.17	Allows registration, recording and viewing of digital video in HD (1280x720) format during EEG registration up to 128 channels. The full version, which contains the ability to view videos with increased and decreased playback

				speeds, and also has a digital zoom.
55	Software option for recording and viewing digital video-EEG up to 256 channels, CLINIC version	SW ACQ DV256	1.17	Allows registration, recording and viewing of digital video in HD format (1280x720) during EEG registration up to 256 channels. The full version, which contains the ability to view videos with increased and decreased playback speeds, and also has a digital zoom.
56	Software option for recording and viewing digital IP video-EEG over a network of up to 128 channels, CLINIC version	SW ACQ DVIP128	1.42	Allows registration, recording and viewing of digital video in HD (1280 x 720) format from cameras connected via IP protocol during EEG registration up to 128 channels. The full version, which contains the ability to view videos with increased and decreased playback speeds, and also has a digital zoom.
57	Software option for recording and viewing digital IP video-EEG over a network of up to 256 channels, CLINIC version	SW ACQ DVIP256	1.42	Allows registration, recording and viewing of digital video in HD (1280 x 720) format from cameras connected via IP protocol during EEG registration up to 256 channels. The full version, which contains the ability to view videos with increased and decreased playback speeds, and also has a digital zoom.

SystemPLUS EVOLUTION - (Test Test (EEG Exam Acquires))


File Modify Acquire Analyse Reference Tools Window Help


0.008 Hz

10 sec


1600 $\mu\text{V}/\text{cm}$

1-All recorded elr

254

SystemPLUS EVOLUTION

Exam Filters: Rec.Date: From 23.04.2021 To

Resource: EEG Patient: Test, Test Exam: EEG Trace

Name	Birthdate	Exam	Date	Stabst	File	Date	Storage	Locked By	Comment
Test, Test	1 Январь 2000 г.	EEG Trace	11.04.2021 4:55:48		EEG Trace	14.04.2021 16:37:55			
		EEG Trace	11.04.2021 12:55:45		EEG Video	14.04.2021 16:37:55			
		EEG Trace	11.04.2021 20:55:42		EEG Video	14.04.2021 17:22:56			
		EEG Trace	12.04.2021 4:55:39		EEG video	14.04.2021 18:07:57			
		EEG Trace	12.04.2021 12:55:37		EEG Video	14.04.2021 18:52:57			
		EEG Trace	12.04.2021 20:55:34		EEG Video	14.04.2021 19:37:58			
		EEG Trace	13.04.2021 4:55:31		EEG Video	14.04.2021 20:22:59			
		EEG Trace	13.04.2021 12:55:28		EEG Video	14.04.2021 21:07:59			
		EEG Trace	13.04.2021 20:55:25		EEG Video	14.04.2021 21:53:00			
		Trace	14.04.2021 11:03:29		EEG Video	14.04.2021 22:38:00			
		Trace	14.04.2021 16:37:55		EEG Video	14.04.2021 23:23:01			
		Trace	15.04.2021 0:37:52		EEG Video	15.04.2021 0:08:02			
		Trace	15.04.2021 8:37:50						
		Trace	15.04.2021 16:37:47						
		Trace	15.04.2021 18:36:39						
		Trace	15.04.2021 18:47:24						

Menu: Main, 10 days, History, Automatic CD/DVD Writing Procedure, View Remote Machine, Remote Machine Take Control, RealTime Networking, Remote Review, RealTime Networking with Control, Delete, Properties

SystemPlus EVOLUTION Working Area Archive

Resource Available GB

Resource	Available GB
Archive2	345,31/2 047,40
EEG	859,04/2 794,50
Local Data	345,31/2 047,40
Secure	345,31/2 047,40
SQLbase	345,31/2 047,40
SQLbase2	345,31/2 047,40
Working Area	345,31/2 047,40
Archive	345,31/2 047,40

Read/Write Local Zone

ENG 14:01 23.04.2021

 ООО ЦРМИ  INFO@CIRMI.RU  WWW.GOSZDRAVNADZOR.RU

Exam: EEG Year: 2021

Exam Filters: Rec. Date: From: To: Mkr: Dr: Th:

F2 F3 F5 F6 F7 F8 F9 F10 F11 F12

EEG EEG EEG EP IOM EEG EEG EEG EEG EEG EEG

Resource: EEG Patient: Test, Test Exam: EEG Trace

Name	Birthdate	Exam	Date	Statist	File	Date	Storage	Locked By	Comment
Test, Test	1 Января 2000 г.	EEG Trace	11.04.2021 4:55:48		EEG Trace	14.04.2021 16:37:55			
		EEG Trace	11.04.2021 12:55:45		EEG Video	14.04.2021 16:37:55			
		EEG Trace	11.04.2021 20:55:42		EEG Video	14.04.2021 17:22:56			
		EEG Trace	12.04.2021 4:55:39		EEG Video	14.04.2021 18:07:57			
		EEG Trace	12.04.2021 12:55:37		EEG Video	14.04.2021 18:52:57			
		EEG Trace	12.04.2021 20:55:34		EEG Video	14.04.2021 19:37:59			
		EEG Trace	13.04.2021 4:55:31		EEG Video	14.04.2021 20:22:59			
		EEG Trace	13.04.2021 12:55:28		EEG Video	14.04.2021 21:07:59			
		EEG Trace	13.04.2021 20:55:25		EEG Video	14.04.2021 21:53:00			
		EEG Trace	14.04.2021 11:03:29		EEG Video	14.04.2021 22:38:00			
		EEG Trace	14.04.2021 16:37:55		EEG Video	14.04.2021 23:23:01			
		EEG Trace	15.04.2021 0:37:52		EEG Video	15.04.2021 0:08:02			
		EEG Trace	15.04.2021 8:37:50						
		EEG Trace	15.04.2021 16:37:47						
		EEG Trace	15.04.2021 18:36:39						
		EEG Trace	15.04.2021 18:47:24						

Properties

Archive2 EEG Local Data secure SQLbase SQLbase2

Working Area Archive

Resource Available GB

Archive2	345,31/2 047,40
EEG	859,04/2 794,50
Local Data	345,31/2 047,40
secure	345,31/2 047,40
SQLbase	345,31/2 047,40
SQLbase2	345,31/2 047,40
Working Area	345,31/2 047,40
Archive	345,31/2 047,40

Read/Write Local Zone

 **ООО ЦПМИ**  **INFO@CPMI.RU**  **WWW.GOSZDRAVNADZOR.RU**

F2

Exam Filters

Rec. Date

From: 04/2021

To:

Resource: EEG

Name	Birthdate
Test, Test	1 Январь 2000 г.

F3 F5 F6 F7 F8

EEG EMG EP IOM INFO

F9 F10 F11 F12

EEG EMG EP RLV

Patient: Test, Test

Exam: EEG Trace

EEG Options

Video Recording Video Reviewing

Headbox EEG EEG Recording Video Source

Trace Length:

☒ No limit
 ☐ Close Trace File Every 1:00
 ☐ Close Trace File At Time 1:00:00
 ☐ Make a Single Trace and Stop Recording

Network Warning Messages

☐ Enable warning text messages

Available Network Resources

Send messages to

NetSend Program

McEvsnReceiver

NetSend.exe

McNetSend.exe

EvmsButtonReceiver

Archive2 EEG Local Data secure

Working Area Archive

SystemPlus EVOLUTION

Date	Storage	Locked By	Command
14.04.2021 16:37:55			
14.04.2021 16:37:55			
14.04.2021 17:22:56			
14.04.2021 18:07:57			
14.04.2021 18:52:57			
14.04.2021 19:37:58			
14.04.2021 20:22:59			
14.04.2021 21:07:59			
14.04.2021 21:53:00			
14.04.2021 22:38:00			
14.04.2021 23:23:01			
15.04.2021 0:08:02			

Available CH	Descr
345,31/2 047,40	
850,04/2 794,50	
345,31/2 047,40	
345,31/2 047,40	
345,31/2 047,40	
345,31/2 047,40	
345,31/2 047,40	
345,31/2 047,40	

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

F2

Exam Filters

Mkr

Th

Rec.Data

From

To

F3 F5 F6 F7 F8

EEG

EMG

EP

IOM

INFO

F9 F10 F11 F12

EEG

EMG

EP

RIP

Resource: EEG

Name

Test, Test

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Patient: Test Test

Exam: EEG Trace

EEG Options

Video Recording

Video Reviewing

Headbox

EEG

EEG Recording

Video Source

Video Acquisition Input Type

Not Present

DV (Firewire, Canopus 1394, Canopus A/DVio)

SD (BlackMagic IntensityPro)

HD (BlackMagic IntensityPro)

SD (BlackMagic IntensityPro 4K)

HD (BlackMagic IntensityPro 4K)

IP Camera (SD)

IP Camera (HD)

Webcam

Settings

Video Capture Source

Microsoft DV Camera and VCR

Select...

Settings

Capture audio

Audio PC Speaker Output

Keep Video Size Dimension

Camera Control

Protocol

Serial Port 1

Video

Palco-D

Sony SNC

Serial Port 2 VAPOR (Axis)

Use external camera control (over IP)

Exit

Help

Archive2

EEG

Local Data

Archive

SystemPlus EVOLUTION

Working Area


Archive

Read/Write

Local Zone

Date	Storage	Locked By	Comment
14.04.2021 16:37:55			
14.04.2021 16:37:55			
14.04.2021 17:22:56			
14.04.2021 18:07:57			
14.04.2021 19:52:57			
14.04.2021 19:37:58			
14.04.2021 20:22:59			
14.04.2021 21:07:59			
14.04.2021 21:53:00			
14.04.2021 22:38:00			
14.04.2021 23:23:01			
15.04.2021 0:08:02			

Available GB	Descr
345,31/2 047,40	
859,04/2 794,50	
345,31/2 047,40	
345,31/2 047,40	
345,31/2 047,40	
345,31/2 047,40	
345,31/2 047,40	
345,31/2 047,40	
345,31/2 047,40	

SystemPlus EVOLUTION

File Edit View Tools Help

F2 Exam Filters Rec.Date From 23.04.2021 To

F3 F5 F6 F7 F8 EEG EMG EP IOM INFO

F9 F10 F11 F12 EEG EMG EP REL

Resource: EEG Patient: Test, Test Exam: EEG Trace

File	Date	Storage	Locked By	Comment
EEG Trace	11.04.2021 4:55:48			
EEG Trace	11.04.2021 12:55:45			
EEG Trace	11.04.2021 20:55:42			
EEG Trace	12.04.2021 4:55:39			
EEG Trace	12.04.2021 12:55:37			
EEG Trace	12.04.2021 20:55:34			
EEG Trace	14.04.2021 16:37:55			
EEG Video	14.04.2021 16:37:55			
EEG Video	14.04.2021 17:22:56			
EEG Video	14.04.2021 18:07:57			
EEG Video	14.04.2021 18:52:57			
EEG Video	14.04.2021 19:37:58			
EEG Video	14.04.2021 20:22:59			
EEG Video	14.04.2021 21:07:59			
EEG Video	14.04.2021 21:53:00			
EEG Video	14.04.2021 22:38:00			
EEG Video	14.04.2021 23:23:01			
EEG Video	15.04.2021 0:08:02			
EEG Trace	15.04.2021 16:37:47			
EEG Trace	15.04.2021 18:36:39			
EEG Trace	15.04.2021 18:47:24			

Message

Export Viewer

☐ Trace

☒ Trace and video

☐ Export selected traces only

OK Cancel

Archive2 EEG Local Data secure SQLbase SQLbase2


Working Area Archive

Resource	Available GB	Descr
Archive2	345,31/2 047,40	
EEG	859,04/2 794,50	
Local Data	345,31/2 047,40	
secure	345,31/2 047,40	
SQLbase	345,31/2 047,40	
SQLbase2	345,31/2 047,40	
Working Area	345,31/2 047,40	
Archive	345,31/2 047,40	1

Read/Write Local Zone

SystemPlus EVOLUTION

Windows taskbar: 16:22 23.04.2021

 Resource: EEG

Name	Birthdate
Test, Test	Saturday, Januar...

👤 Patient: Test, Test

Exam

 EEG Trace

FFC T-260

100

ILG 220 11308

EEG Trace

 EEG Trace

EEG Trace

200

EEC

FIG. 1 EEG trace

EEG Trace

EEG Trace

Call 855-333-3333

E.E.C.

 EEG Trace EEG Trace

EEG Trace

110

146

EEG Trace

 EEG Trace

EEG Trace

The screenshot displays the NeuroLab software interface with several settings panels:

- Video Recording:** Includes a play button icon and the text "Video Recording".
- Headbox:** Includes an EEG icon and the text "EEG Recording".
- Video Sources:** Includes a camera icon and the text "Video Sources".
- Video Acquisition Input Type:**
 - ☒ Not Present
 - ☐ SD (BlackMagic IntensityPro)
 - ☐ SD (BlackMagic IntensityPro) (HDMI)
 - ☐ HD (BlackMagic IntensityPro)
 - ☐ HD (BlackMagic IntensityPro) (HDMI)
 - ☐ IP Camera (SD)
 - ☐ IP Camera (HD)
 - ☐ Webcam
- Video Capture Source:**
 - ☒ Decklink Video Capture
 - ☐ Select...
- Capture audio:**
 - ☒ Decklink Audio Capture
 - ☐ Audio PC Speaker Output
- Keep Video Size Dimension:** ☐
- Camera Control:**
 - Protocol:
 - Serial Port 1:
 - Serial Port 2:
 - ☐ Use external camera control (over IP)
- Video Preview:** A small window showing a live video feed of a person's face.

Exam: EEG Trace

	Date	Storage	Locked By	Comment
EEG Trace	4/23/2021 10:46:22 ..			

The screenshot shows the 'Properties' dialog box for 'Video Format'. The 'Video Format' dropdown is set to '1080i50 8 bit 4:2:2 YUV'. The 'Timecode Format' dropdown is set to 'Timecode'. Under 'Timecode Format', the 'Read from VITC (SHIPT L28)' radio button is selected, and the 'Read from MANC (SHIPT RP 188)' radio button is unselected. The 'Blackmagic design' logo is visible in the background.

SystemPlus EVOLUTION

Resource


EEG
Local Data sq
Working Area

849.92/2,794.50
345.27/2,047.40

Descr

SystemPLUS EVOLUTION

File Modify Record Analyse Reference Tools Move Window Help

150

Measure Cursor

EEGEvents Setup

Maps

Data Average

EEG Analyser

IMR1 Artifact Suppression

Analysis Manager

Reveal Detector

Reveal event notification

Smart EEG

Export wide band EEG

Persyst Enable

Persyst Settings

Manual hypnogram

Crisis

External Program

External Program Setup

Digital Trigger Insertion Parameters

Analog Trigger Insertion Parameters

Display trigger information

Trigger deletion

Insert trigger on muscular activity

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Timer: 03:56:21.0
Hour: 14:42:43
Date: 23.04.2021


0.008 Hz
10 sec
200 μ V/cm
1-All recorded data
Zoom x 1.0 Speed x 1.0

14:42:43.00

Volume

Speed

Patient: Test_Test

Exam: EEG Trace

Name	Birthdate
Test, Test	Saturday, Januar...

[illegible]

Video Recording

EEG Recording

Video Source

Video Acquisition Input Type

Not Present

SD (BlackMagic IntensityPro)

HD (BlackMagic IntensityPro)

SD (BlackMagic IntensityPro 4K1G)

HD (BlackMagic IntensityPro 4K1G)

IP Camera (SD)

IP Camera (HD)

Webcam

Video Capture Source

Decklink Video Capture

Capture audio

Decklink Audio Capture

Audio PC Speaker Output

Keep Video Size Dimension

Camera Control

Protocol

Serial Port 1

COM2

Serial Port 2

(none)

Use external camera control (over IP)

Properties

Audio Format | Audio Stream

Audio Format

2 channels, 48kHz, 16-bit

	Date	Storage	Locked By	Comment
EEG Trace	4/23/2021 10:46:22 ...			
EEG Video	4/23/2021 10:46:22 ...			
EEG Video	4/23/2021 11:31:23 ...			
EEG Video	4/23/2021 12:16:24 ...			
EEG Video	4/23/2021 1:01:24 PM			
EEG Video	4/23/2021 1:46:25 PM			
EEG Video	4/23/2021 2:31:25 PM			
EEG Video	4/23/2021 3:16:26 PM			
EEG Video	4/23/2021 4:01:27 PM			
EEG Video	4/23/2021 4:46:27 PM			
EEG Video	4/23/2021 5:31:28 PM			

Resource	
EEG	
Local Data 501	
Working Area	

Properties

Audio Format | Audio Stream

Blackmagic design



Audio Format

2 channels, 48KHz, 16-bit

OK Cancel Apply

SystemPLUS EVOLUTION - (Test Test (EEG Exam Acquire))

File Modify Acquire Analyse Reference Tools Move Window Help

0.008 Hz


10 sec

1600 $\mu\text{V}/\text{cm}$


1-All recorded ch

Test, Test Monday, April 26, 2021 12:15 (View Eeg Trace)

Test, Test (Eeg Export Acquire)

Measure Cursor

EEGEvents Setup

Maps

Amplitude Map

Data Average

Frequency Maps

EEG Analyser

Coherence Maps

fMRI Artifact Suppression

Analysis Manager

Reveal Detector

Reveal event notification

Smart EEG

Export wide band EEG

Peryst Enable

Peryst Settings

Manual hypnogram

Critique

External Program

External Program Setup

Fp1 G2

F7 G2

F3 G2

F1 G2

F4 G2

F8 G2

T3 G2

C3 G2

C2 G2

C4 G2

T4 G2

T5 G2

P3 G2

Pz G2

P4 G2

T6 G2

O1 G2

O2 G2

eA22

eA23

eA24

eA25

eA26

eA27

eA28

eA29

eA30

eA31

eA32

eB01

eB02

eB03

eB04

eB05

eB06

eB07

eB08

eB09

eB10

eB11

eB12

eB13

eB14

eB15

eB16

eB17

eB18

eB19

eB20

eB21

eB22

eB23

eB24

eB25

eB26

eB27

eB28

eB29

eB30

eB31

eB32

MGR - MGR

Timer: 00:00:00

Hour: 10:47:01

Date: 18.04.2021

SystemPLUS

1,600 Hz 70 Hz 15 sec 100 μ V/cm 1-Reduced

fMRI Artifact Suppression

Gradient Artifact

fMRI Acquisition Parameters

Number of Total Volumes: 200 TR (msec): 3000

Exclude Volumes from Average:

Filtering Algorithm

☒ Averaging only ☐ LMS Adaptive Algorithm ☐ NLMS Adaptive Algorithm ☐ RLS Adaptive Algorithm

Algorithm Parameters

Reference Volumes: 12 Step Size Divisor (LMS): 20

Adaptive Filter Order: 30 μ (NLMS): 0.15

Lambda (RLS): 0.98

☒ Volume Realignment ☒ Two Step Filtering

Auto-align Volumes

Numbering Volumes

Artifact Filtering

Cardioballogram Artifact


Cardio Electrode Position: PC Peak Delay (msec): 200

Identify BRS

Polynomial Subtraction

Exit

C64 - 258Hz

Select

Video Enable

View Impedance

Filter Trace

Archive Montage

Options

Duplicate and lock

Acquisition Video Split

Acquisition Video Split Refresh

As Recorded


Stimulator Control

Stimulator Program Setup

LTM Cortical Stim

LTM Cortical Stim Setup

EEG Sound

Timer: 00:00:00.0


Hour : 12:27:51


Date : 4/26/2021

Type here to search

12:31 PM 4/26/2021 ENG

0.008 Hz


10 sec

200 μ V/cm

1-All recorded el

Select


Display Base Time[00:40:00] Adapt Time Base to Trace: Off

Type here to search

12:54 PM 4/26/2021 ENG

Exam Filters

From: 15/02/2018 To:

Rec Date:

F2 F3 F4 F5 F6 F7 F8 F9 F10 F11 F12

INFO

REP

Resource: BrainQuick

Name	Birthday	Address	City	Province	Exam	Date	Statistic	File	Date	Comment	Storage	Locked By
------	----------	---------	------	----------	------	------	-----------	------	------	---------	---------	-----------

Savina, Delta	03/04/2004				EEG	08/05/2018 13:08:45		EEG Trace	08/05/2018 13:08:45			BrainQuick
Criciuth, Tommaso	01/11/2005				EEG			EEG Trace				BrainQuick
Romanini, Martina	10/07/2006							REP New View				BrainQuick
Cardello, Matteo	11/03/2010							EEG Copy To...				BrainQuick
Vittuzzi, Maria R...	01/01/1944							Reduce				
Garbagnoli, Renato	22/11/1935							Archived Montages Propagation				
Sani, Gaetano	03/06/2000							Statistic Items				
Micromed, Feg	23/01/1993							Analyse				
Test, Ee	04/04/1992							Export Viewer				
Pesaresi, Valeria	28/12/1998							Modify				
Petkov, P.	24/12/1940							X Remove				
Biondi, H	21/12/1941							Properties				
Frequenza, Diverse	24/02/2017											
Sab, Alexander	27/12/1991											
Osas, Compilata	01/05/1947											
P. Matia	08/06/1990											
Zeh, Patient	01/01/1900											
User Cld, Patient	01/01/1900											
Case 2, Samuel	15/11/1989											
Intensity Pro, 4k	10/04/2018											
High Def, Video	23/09/2007											
Scotti, Filippo G...	14/04/2002											
Test, Cx000	13/03/2017											
Gay, Bernard	02/10/1939											
Casati, Emily	09/05/1991											
Ful Hd, John	27/05/2003											
Cfm, Milano	25/08/2006											
Rozzo, Nicola	25/04/1971											
Ferrari, Carlo	11/03/1955											

Patient: Scotti, Filippo Giulio

Exam: EEG Trace

File	Date	Comment	Storage	Locked By
------	------	---------	---------	-----------

EEG	08/05/2018 13:08:45			BrainQuick
REP				BrainQuick
EEG Copy To...	08:45			BrainQuick
Reduce				
Archived Montages Propagation				
Statistic Items				
Analyse				
Export Viewer				
Modify				
X Remove				
Properties				

micromed

BrainQuick

Cardinal

ECG/EMG

HD EXT

2m

EEG view

Viewer

Video

VIEWING

VERSION

Resource

BrainQuick

Cardinal

DEMO EKG

HD EXT

Available on

2022/11/06/14


Not available

2022/11/06/14

Not available


2022/11/06/14

HIS Workflow Management with Plug-In

micromed


HIS Workflow Management with Plug-In




micromed

165

HIS Workflow Management with Plug-In

micromed

SystemPLUS EVOLUTION

F2

Report

Options...

Ambulatory recorder

EMG Quick Key

Open New Report Model

Backlog Schedule

Delete/Modify Archive Data

Rec.Date

from 26.04.2021

Memory card reader

Direct acquisition

Setup

F3 F5 F6 F7 F8

EEG EMG EP IOM INFO

F9 F10 F11 F12

EEG EMG EP REP

Patient: Test, Test

Exam: EEG Trace

Exam	Date	Statist	File	Date	Storage	Locked By	Comment
EEG Trace	18.04.2021 10:47:01		EEG Trace	23.04.2021 10:46:22			
EEG Trace	18.04.2021 18:46:58		EEG Video	23.04.2021 10:46:22			
EEG Trace	19.04.2021 2:46:56		EEG Video	23.04.2021 11:31:23			
EEG Trace	19.04.2021 10:46:53		EEG Video	23.04.2021 12:16:24			
EEG Trace	19.04.2021 18:46:50		EEG Video	23.04.2021 13:01:24			
EEG Trace	20.04.2021 2:46:46		EEG Video	23.04.2021 13:46:25			
EEG Trace	20.04.2021 10:46:44		EEG Video	23.04.2021 14:31:25			
EEG Trace	20.04.2021 18:46:41		EEG Video	23.04.2021 15:16:26			
EEG Trace	21.04.2021 2:46:39		EEG Video	23.04.2021 16:01:27			
EEG Trace	21.04.2021 10:46:38		EEG Video	23.04.2021 16:46:27			
EEG Trace	21.04.2021 18:46:36		EEG Video	23.04.2021 17:31:28			
EEG Trace	22.04.2021 2:46:33						
EEG Trace	22.04.2021 10:46:30						
EEG Trace	22.04.2021 18:46:28						
EEG Trace	23.04.2021 2:46:25						
EEG Trace	23.04.2021 10:46:22						

Archive2

EEG

Local Data

secure

SQLbase

SQLbase2

Working Area

Archive

Resource

Available GB

Descr

Archive2

EEG

Local Data

secure

SQLbase

SQLbase2

Working Area

Archive

345,26/2 047,40

345,26/2 047,40

345,26/2 047,40

345,26/2 047,40

345,26/2 047,40

345,26/2 047,40

345,26/2 047,40


Read/Write

Local Zone


ENG

26.04.2021

269

Select

Timer 00:00:00
Hour 10:46:22
Date 23.04.2021

270

SystemPLUS EVOLUTION

Exam Filters: Mkr Dr Th. Rec.Date: From 26.04.2021 To


Resource: EEG Patient: Test, Test Exam: EEG Trace


Name	Birthdate	Exam	Date	Stabst	File	Date	Storage	Locked By	Comment
Test, Test	1.04.2021	EEG Trace	18.04.2021 10:47:01		EEG Trace	23.04.2021 16:01:27			
		EEG Trace	18.04.2021 18:46:58		EEG Trace	23.04.2021 16:46:27			
		EEG Trace	19.04.2021 2:46:56		EEG Trace	23.04.2021 17:31:28			
		EEG Trace	19.04.2021 10:46:53						
		EEG Trace	19.04.2021 18:46:50						
		EEG Trace	20.04.2021 2:46:46						
		EEG Trace	20.04.2021 10:46:44						
		EEG Trace	20.04.2021 18:46:41						
		EEG Trace	21.04.2021 2:46:39						
		EEG Trace	21.04.2021 10:46:38						
		EEG Trace	21.04.2021 18:46:36						
		EEG Trace	22.04.2021 2:46:33						
		EEG Trace	22.04.2021 10:46:30						
		EEG Trace	22.04.2021 18:46:28						
		EEG Trace	23.04.2021 2:46:25						
		EEG Trace	23.04.2021 10:46:22						

SystemPlus EVOLUTION Working Area Archive

Resource	Available GB	Descr
Archive2	345,26/2 047,40	
EEG	849,89/2 794,50	
Local Data	345,26/2 047,40	
Secure	345,26/2 047,40	
SQLbase	345,26/2 047,40	
SQLbase2	345,26/2 047,40	
Working Area	345,26/2 047,40	
Archive	345,26/2 047,40	1

Read/Write Local Zone

F2

Exam Filters

From

To

Rec Date

26.04.2021

Resource: EEG

Name	Birthdate
Test, Test	1 Январь 2000 г.

F3

F5

F6

F7

F8

EEG

EMG

EP

IOM

INFO

F9

F10

F11

F12

EEG

EMG

EP

REP

Exam: EEG Trace

File	Date	Storage	Locked By	Comment
EEG Trace	23.04.2021 10:46:22			
EEG Video	23.04.2021 10:46:22			
EEG Video	23.04.2021 11:31:23			
EEG Video	23.04.2021 12:16:24			
EEG Video	23.04.2021 13:01:24			
EEG Video	23.04.2021 13:46:25			
EEG Video	23.04.2021 14:31:25			
EEG Video	23.04.2021 15:16:26			
EEG Video	23.04.2021 16:01:27			
EEG Video	23.04.2021 16:46:27			
EEG Video	23.04.2021 17:31:28			

Available GB

1042

23.04.2021

Video Recording

Video Reviewing

Headbox

EEG Recording

Video Source

Video Acquisition Input Type

Not Present

SD (BlackMagic IntensityPro)

SD (BlackMagic IntensityPro 4K)

IP Camera (SD)

Webcam

DV (Firewire, Canopus 1394, Canopus AEDVio)

HD (BlackMagic IntensityPro)

HD (BlackMagic IntensityPro 4K)

IP Camera (HD)

Sony SNC-EP521, EP580, EP581

Video Capture Source

IP Camera Address

Capture audio

Audio PC Speaker Output

Keep Video Size Dimension

Discover

Settings

Select...

Settings

Camera Control

Protocol

Serial Port 1

Serial Port 2

Use external camera control (over IP)

Right in the Intendence Check

Exit

Help

Archive2

EEG

Local Data

secure

SystemPlus EVOLUTION

Working Area

Archive

Read/Write

Local Zone

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

Resource: EEG	Patient: Test, Test	Exam: EEG Trace
---------------	---------------------	-----------------

 Archive2
  EEG
  Local Data
  secure
  SQLbase
  SQLbase2

 Working Area
  Archive

Resource	Available GB	Descr
Archive2	345,26/2 047,40	
EEG	849,89/2 794,50	
Local Data	345,26/2 047,40	
secure	345,26/2 047,40	
SQLbase	345,26/2 047,40	
SQLbase2	345,26/2 047,40	
Working Area	345,26/2 047,40	
Archive	345,26/2 047,40	

SystemPLUS EVOLUTION - [Test Test (EEG Exam Acquire)]

Modify Acquire Analyse Reference Tools Move Window Help

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
High pass filter: 0.008 Hz Low pass filter: 1000 Hz Gain: 1600 $\mu\text{V}/\text{cm}$



Reading data from cable


C64 - 258Hz


Power 81%

Type here to search

Block diagram of connecting the components of the system. Assembled.

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AUTHENTICATION OF SIGNATURE

I the undersigned Nicoletta Mammucari, Notary public in Mogliano Veneto, duly enrolled in the Notarial College of the District of Treviso, who knows English language, attest the authenticity of the signature of Mr. Cristiano RIZZO, born in Venice on the 2th of March 1969, resident for the office as follows, as special attorney of the Company **MICROMED SPA**, with registered office in Mogliano Veneto (TV) Via Giotto n.2, tax code and registration number in the Treviso - Belluno Company Register n.03906850262, right power of attorney of 24th of July 2018, duly filed in the company register, whose personal identity, qualification and powers I, the Notary, am certain.
Mogliano Veneto, Via IV Novembre n.6/1, on the fourteenth of February twothousandtwentythree.

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



145

SystemPLUS EVOLUTION

Exam Filters: [Patient Name] Rec.Date: [Date]

Resource: EEG Patient: Test, Test Exam: EEG Trace

Name	Birthdate	Exam	Date	Statist	File	Date	Storage	Locked By	Comment
Test, Test	1 Января 2000 г.	EEG Trace	11.04.2021 4:55:48		EEG Trace	14.04.2021 16:37:55			
		EEG Trace	11.04.2021 12:55:45		EEG Video	14.04.2021 16:37:55			
		EEG Trace	11.04.2021 20:55:42		EEG Video	14.04.2021 17:22:56			
		EEG Trace	12.04.2021 4:55:39		EEG Video	14.04.2021 18:07:57			
		EEG Trace	12.04.2021 12:55:37		EEG Video	14.04.2021 18:52:57			
		EEG Trace	12.04.2021 20:55:34		EEG Video	14.04.2021 19:37:58			
		EEG Trace	13.04.2021 4:55:31		EEG Video	14.04.2021 20:22:59			
		EEG Trace	13.04.2021 12:55:28		EEG Video	14.04.2021 21:07:59			
		EEG Trace	13.04.2021 20:55:25		EEG Video	14.04.2021 21:53:00			
		Trace	14.04.2021 11:03:29		EEG Video	14.04.2021 22:38:00			
		Trace	14.04.2021 16:37:55		EEG Video	14.04.2021 23:23:01			
		Trace	15.04.2021 0:37:52		EEG Video	15.04.2021 0:08:02			
		Trace	15.04.2021 8:37:50						
		Trace	15.04.2021 16:37:47						
		Trace	15.04.2021 18:36:39						
		Trace	15.04.2021 18:47:24						

Properties

Archive2 EEG Local Data secure SQLbase SQLbase2

Working Area Archive

Resource	Available GB
Archive2	145,31/2 047,40
EEG	859,04/2 794,50
Local Data	345,31/2 047,40
secure	345,31/2 047,40
SQLbase	345,31/2 047,40
SQLbase2	345,31/2 047,40
Working Area	345,31/2 047,40
Archive	345,31/2 047,40

SystemPlus EVOLUTION

Windows Taskbar: 16:03 23.04.2021

SystemPLUS EVOLUTION

Exam Filters Rec.Date

Resource: EEG Patient: Test, Test Exam: EEG Trace

Name	Birthdate	Exam	Date	Statist	File	Date	Storage	Locked By	Comment
Test, Test	1 Январь 2000 г.	EEG Trace	11.04.2021 4:55:48		EEG Trace	14.04.2021 16:37:55			
		EEG Trace	11.04.2021 12:55:45		EEG Video	14.04.2021 16:37:55			
		EEG Trace	11.04.2021 20:55:42		EEG Video	14.04.2021 17:22:56			
		EEG Trace	12.04.2021 4:55:39		EEG Video	14.04.2021 18:07:57			
		EEG Trace	12.04.2021 12:55:37		EEG Video	14.04.2021 18:52:57			
		EEG Trace	12.04.2021 20:55:34		EEG Video	14.04.2021 19:37:58			
		EEG Trace	13.04.2021 4:55:31		EEG Video	14.04.2021 20:22:59			
		EEG Trace	13.04.2021 12:55:28		EEG Video	14.04.2021 21:07:59			
		EEG Trace	13.04.2021 20:55:25		EEG Video	14.04.2021 21:53:00			
		Trace	14.04.2021 11:03:29		EEG Video	14.04.2021 22:38:00			
		Trace	14.04.2021 16:37:55		EEG Video	14.04.2021 23:23:01			
		Trace	15.04.2021 0:37:52		EEG Video	15.04.2021 0:08:02			
		Trace	15.04.2021 8:37:50						
		Trace	15.04.2021 16:37:47						
		Trace	15.04.2021 18:36:39						
		Trace	15.04.2021 18:47:24						

Archive2 EEG Local Data secure SQLbase SQLbase2

Working Area Archive

SystemPlus EVOLUTION

Read/Write Local Zone

Resource Available GB Descr

Archive2	345,31/2 047,40	
EEG	859,04/2 794,50	
Local Data	345,31/2 047,40	
secure	345,31/2 047,40	
SQLbase	345,31/2 047,40	
SQLbase2	345,31/2 047,40	
Working Area	345,31/2 047,40	
Archive	345,31/2 047,40	

SystemPLUS EVOLUTION

Exam Filters Rec Date From To

Resource: EEG Patient: Test, Test Exam: EEG Trace

File	Date	Storage	Locked By	Comment
EEG Trace	11.04.2021 4:55:48			
EEG Trace	11.04.2021 12:55:45			
EEG Trace	11.04.2021 20:55:42			
EEG Trace	12.04.2021 4:55:39			
EEG Trace	12.04.2021 12:55:37			
EEG Trace	12.04.2021 20:55:34			
EEG Trace	13.04.2021 4:55:31			
EEG Trace	13.04.2021 12:55:28			
EEG Trace	13.04.2021 20:55:25			
EEG Trace	14.04.2021 11:03:29			
EEG Trace	14.04.2021 16:37:55			
EEG Trace	15.04.2021 0:37:52			
EEG Trace	15.04.2021 8:37:50			
EEG Trace	15.04.2021 16:37:47			
EEG Trace	15.04.2021 18:36:39			
EEG Trace	15.04.2021 18:47:24			
EEG Trace	14.04.2021 16:37:55			
EEG Video	14.04.2021 16:37:55			
EEG Video	14.04.2021 17:22:56			
EEG Video	14.04.2021 18:07:57			
EEG Video	14.04.2021 18:52:57			
EEG Video	14.04.2021 19:37:58			
EEG Video	14.04.2021 20:22:59			
EEG Video	14.04.2021 21:07:59			
EEG Video	14.04.2021 21:53:00			
EEG Video	14.04.2021 22:38:00			
EEG Video	14.04.2021 23:23:01			
EEG Video	15.04.2021 0:08:02			

Properties

Archive2 Local Data secure SQLbase SQLbase2

Working Area Archive

Resource	Available GB	Descr
Archive2	345,31/2 047,40	
EEG	859,04/2 794,50	
Local Data	345,31/2 047,40	
secure	345,31/2 047,40	
SQLbase	345,31/2 047,40	
SQLbase2	345,31/2 047,40	
Working Area	345,31/2 047,40	
Archive	345,31/2 047,40	

Read/Write Local Zone

Resource: EEG	Patient: Test, Test	Exam: EEG Trace
---------------	---------------------	-----------------

[illegible]

The Role of the Teacher

- The teacher is a facilitator of learning
- The teacher is a guide and supporter
- The teacher is a resource
- The teacher is a learner
- The teacher is a collaborator
- The teacher is a reflective practitioner
- The teacher is a professional
- The teacher is a leader
- The teacher is a change agent
- The teacher is a lifelong learner

The screenshot shows a presentation window titled "Presentation" with a standard toolbar at the top. The slide content includes a title "The Role of the Teacher" and a bulleted list of ten roles. The roles listed are: The teacher is a facilitator of learning, The teacher is a guide and supporter, The teacher is a resource, The teacher is a learner, The teacher is a collaborator, The teacher is a reflective practitioner, The teacher is a professional, The teacher is a leader, The teacher is a change agent, and The teacher is a lifelong learner.

SystemPlus Evolution

Exam Filters

From 23.04.2021 To

Rec. Date

Source: EEG

Birthdate

19.04.2000 F.

Test

Patient: Test, Test

EEG Options

Video Recording

Headbox EEG

Trace Length

No limit

Close Trace File Every:

Close Trace File At Time:

Make a Single Trace and Stop Recording

Network Warning Messages

Enable warning text messages

Available Network Resources

Send messages to

Refresh

NetSend Program

McEventReceiver

NetSend.exe

McNetSend.exe

MS-Exchange-2003

EventSubscriptionReceiver

Exam: EEG Trace

Storage

Locked By

Date

14.04.2021 16:37:55
14.04.2021 16:37:55
14.04.2021 17:22:56
14.04.2021 18:07:57
14.04.2021 18:52:57
14.04.2021 19:37:58
14.04.2021 20:22:59
14.04.2021 21:07:59
14.04.2021 21:53:00
14.04.2021 22:38:00
14.04.2021 23:23:01
15.04.2021 0:08:02

Archiv2

EEG

Local Data

secure

Working Area

Archive

SystemPlus EVOLUTION

Read/Write


Local Zone

Available GB	Descr
345,31/2 047,40	1
859,04/2 794,50	
345,31/2 047,40	
345,31/2 047,40	
345,31/2 047,40	
345,31/2 047,40	
345,31/2 047,40	

16.11

15.04.2021

546

SystemPlus EVOLUTION

Exam Filters: Rec.Date: 23.04.2021

F2 F3 F5 F6 F7 F8 F9 F10 F11 F12

EEG EMG EP IOM INFO EEG EMG EP REP

Resource: EEG Patient: Test, Test Exam: EEG Trace

Name	Birthdate	Exam	Date	Statist	File	Date	Storage	Locked By	Comment
Test, Test	1 Январь 2000 г.	EEG EEG Trace	11.04.2021 4:55:48		EEG EEG Trace	14.04.2021 16:37:55			

General Settings Jobs Collection

TransferJob Edit

Job Name: Transfer Job Description: Enabled

Source host domain credentials

Host name: Domain: WORKGROUP User name: Password: Connect to Host and read Resources

Local resource database folder: c:\system98\settings

Transfer type: Copy to Security and move to Working Area

Source Resource Copy Settings Move Settings

Resource name: SQLbase2 Resource Path: Browse...

Data source: Microsoft SQL Server (SQL Client) Server name: Log on to the server: Use Windows Authentication Use SQL Server Authentication

User name: Password: Database name: base Test connection

Validate Cancel OK

Working Area Archive 345,31/2 047,40 345,31/2 047,40 1

Read/Write Local Zone

File Edit View Tools Help

F2 F3 F5 F6 F7 F8 F9 F10 F11 F12

Exam Filters Rec.Date

Resource: EEG Patient: Test, Test Exam: EEG Trace

Name	Birthdate	Exam	Date	Storage	Locked By	Comment
Test, Test	Saturday, Januar...	EEG Trace	4/23/2021 10:46:22 ...			

EEG Trace

Video Recording Video Reviewing

Headbox EEG Recording Video Source

Video Acquisition Input Type

Not Present DV (Firewire, Camopus 1394, Camopus ADE DVio)

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

Headbox EEG EEG Recording Video Source

Video Acquisition Input Type

☒ Not Present ☐ DV (Firewire, Camcorder 1394, Camcorder ACE DVio)

☐ SD (BlackMagic IntensityPro) ☐ HD (BlackMagic IntensityPro)

☐ SD (BlackMagic IntensityPro 4K) ☐ HD (BlackMagic IntensityPro 4K)

☐ IP Camera (SD) ☐ IP Camera (HD)

☐ Webcam

Video Capture Source

☐ Decklink Video Capture

☒ Capture audio

☐ Decklink Audio Capture

☐ Audio PC Speaker Output

☐ Keep Video Size Dimension

Camera Control

Protocol

☐ Use external camera control (over IP)

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





Adaptive Filter Order: 20 LMS (RLMS): 0.15

Lambda (RLS): 0.5

☒ Volume reassignment ☒ Fast Noise Removal

Cardioballogram Artifact

Cardio Electrode Position: Pz Peak Delay (msec): 200

Fp2 C4

C4 O2

Fp2 T4

T4 O2

Fp1 C3


C3 O1

Fp1 T3


T3 O1

MKR+ MKR-


C64 - 256Hz


Select


- Video Enable
- View Impedance
- Filter Trace
- Archive Montage
- Options
- Duplicate and lock
- Acquisition Video Split
- Acquisition Video Split Refresh
- As Recorded
- Stimulator Control
- Stimulator Program Setup
- LTM Cortical Stim
- LTM Cortical Stim Setup
- EEG Sound

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






micromed

SystemPLUS EVOLUTION - [Test Test Friday, April 23, 2021 10:46 (View Log Trace)]

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File Edit View Tools Help

F2 F3 F5 F6 F7 F8 F9 F10 F11 F12

Import Options... Ambulatory recorder EMG Quick Key Open New Report Model Backup Scheduler Statistic Amount Delete/Modify Archive Data

Rec.Date 26.04.2021

Memory card reader Direct acquisition Setup

EEG EMG EP IOM INFO

Patient: Test, Test

Exam: EEG Trace

Name	Exam	Date	Statist	File	Date	Storage	Locked By	Comment
Test, Test	EEG Trace	18.04.2021 10:47:01		EEG Trace	23.04.2021 10:46:22			
	EEG Trace	18.04.2021 10:46:58		EEG Video	23.04.2021 10:46:22			
	EEG Trace	19.04.2021 2:46:56		EEG Video	23.04.2021 11:31:23			
	EEG Trace	19.04.2021 10:46:53		EEG Video	23.04.2021 12:16:24			
	EEG Trace	19.04.2021 18:46:50		EEG Video	23.04.2021 13:01:24			
	EEG Trace	20.04.2021 2:46:46		EEG Video	23.04.2021 13:46:25			
	EEG Trace	20.04.2021 10:46:44		EEG Video	23.04.2021 14:31:25			
	EEG Trace	20.04.2021 18:46:41		EEG Video	23.04.2021 15:16:26			
	EEG Trace	21.04.2021 2:46:39		EEG Video	23.04.2021 16:01:27			
	EEG Trace	21.04.2021 10:46:38		EEG Video	23.04.2021 16:46:27			
	EEG Trace	21.04.2021 18:46:36		EEG Video	23.04.2021 17:31:28			
	EEG Trace	22.04.2021 2:46:33						
	EEG Trace	22.04.2021 10:46:30						
	EEG Trace	22.04.2021 18:46:28						
	EEG Trace	23.04.2021 2:46:25						
	EEG Trace	23.04.2021 10:46:22						

SystemPlus EVOLUTION

Archive2 EEG Local Data secure SQLbase SQLbase2


Working Area Archive

Resource	Available GB	Descr
Archive2	345,26/2 047,40	
EEG	849,89/2 794,50	
Local Data	345,26/2 047,40	
secure	345,26/2 047,40	
SQLbase	345,26/2 047,40	
SQLbase2	345,26/2 047,40	
Working Area	345,26/2 047,40	
Archive	345,26/2 047,40	

Read/Write Local Zone

SystemPLUS EVOLUTION - [Test Test 23 Апрель 2021 10:46 (View Eeg Trace)]

File Modify Analyse Reference Tools Menu Window Help

1.600 Hz 70 Hz 15 sec 100 μ V/cm 1-Reduced

35

B / U T

1 2 3 4 5 6 7 8 9 10

Select

Exam Edit View Tools Help

F2 F3 F5 F6 F7 F8 F9 F10 F11 F12
 EEG EEG EEG EEG EEG EEG EEG EEG EEG EEG EEG EEG
 EEG EEG EEG EEG EEG EEG EEG EEG EEG EEG EEG EEG

Exam Filters Rec.Date
 Mkr: From: 26.04.2021 To:

Resource: EEG Patient: Test, Test Exam: EEG Trace

Name	Birthdate	Exam	Date	Statist	File	Date	Storage	Locked By	Comment
Test, Test	1 Января 2000 г.	EEG EEG Trace	18.04.2021 10:47:01		EEG	23.04.2021 16:15:27			
		EEG EEG Trace	18.04.2021 18:46:58		EEG				
		EEG EEG Trace	19.04.2021 2:46:56		EEG				
		EEG EEG Trace	19.04.2021 10:46:53		EEG				
		EEG EEG Trace	19.04.2021 18:46:50		EEG				
		EEG EEG Trace	20.04.2021 2:46:46		EEG				
		EEG EEG Trace	20.04.2021 10:46:44		EEG				
		EEG EEG Trace	20.04.2021 18:46:41		EEG				
		EEG EEG Trace	21.04.2021 2:46:39		EEG				
		EEG EEG Trace	21.04.2021 10:46:38		EEG				
		EEG EEG Trace	21.04.2021 18:46:36		EEG				
		EEG EEG Trace	22.04.2021 2:46:33		EEG				
		EEG EEG Trace	22.04.2021 10:46:30		EEG				
		EEG EEG Trace	22.04.2021 18:46:28		EEG				
		EEG EEG Trace	23.04.2021 2:46:25		EEG				
		EEG EEG Trace	23.04.2021 10:46:22		EEG				


EEG EEG Video 23.04.2021 16:01:27
 EEG EEG Video 23.04.2021 16:46:27
 EEG EEG Video 23.04.2021 17:31:28


Archive2 EEG Local Data secure SQLbase SQLbase2
 Working Area Archive

SystemPlus EVOLUTION

Resource Available GB Descr
 Archive2 345,26/2 047,40
 EEG 849,89/2 794,50
 Local Data 345,26/2 047,40
 secure 345,26/2 047,40
 SQLbase 345,26/2 047,40
 SQLbase2 345,26/2 047,40
 Working Area 345,26/2 047,40
 Archive 345,26/2 047,40

Read/Write Local Zone

File View Tools Help

F2 F3 F5 F6 F7 F8 F9 F10 F11 F12

Exam Filters Rec.Date

From 26.04.2021 To

Resource: EEG

Patient: Test, Test

Exam: EEG Trace

EEG Options

Video Recording Video Reviewing

Headbox EEG Recording Video Source

Video Acquisition Input 1 type

☒ Not Present ☐ DV (Firewire, Canopus 1394, Canopus A/D EVO)

☐ SD (BlackMagic IntensityPro) ☐ HD (BlackMagic IntensityPro)

☐ SD (BlackMagic IntensityPro 4K) ☐ HD (BlackMagic IntensityPro 4K)

☒ IP Camera (ISD) ☐ IP Camera (MD)

☐ Webcam

Video Capture Source

☒ IP Camera Address: Discover Settings

☐ Capture audio Select Settings

☐ Audio PC Speaker Output

☐ Keep Video Size Dimension

Camera Control

Protocol: Disabled

Send Port 1: none

Send Port 2: none

☐ Use external camera control (over IP)

Exit Help

Working Area Archive

Read/Write Local Zone

File	Date	Storage	Locked By	Comment
EEG Trace	23.04.2021 10:46:22			
EEG Video	23.04.2021 10:46:22			
EEG Video	23.04.2021 11:31:23			
EEG Video	23.04.2021 12:16:24			
EEG Video	23.04.2021 13:31:24			
EEG Video	23.04.2021 13:46:25			
EEG Video	23.04.2021 14:31:25			
EEG Video	23.04.2021 15:16:26			
EEG Video	23.04.2021 16:01:27			
EEG Video	23.04.2021 16:46:27			
EEG Video	23.04.2021 17:31:28			

Available GB Descr

345,26/2 047,40

840,89/2 794,50

345,26/2 047,40

345,26/2 047,40

345,26/2 047,40

345,26/2 047,40

345,26/2 047,40

345,26/2 047,40

345,26/2 047,40

SystemPlus EVOLUTION

Archive2 EEG Local Data secure

Working Area Archive

SystemPLUS EVOLUTION interface showing patient data and EEG traces.

Exam Filters | **Rec.Data**



Resource: EEG | **Patient: Test, Test** | **Exam: EEG Trace**

Name	Birthdate	Exam	Date	Statist	File	Date	Storage	Locked By	Comment
Test, Test	1 Января 2000 г	EEG Trace	18.04.2021 10:47:01		EEG Trace	23.04.2021 10:46:22			
		EEG Trace	18.04.2021 18:46:58		EEG Video	23.04.2021 10:46:22			
		EEG Trace	19.04.2021 2:46:56		EEG Video	23.04.2021 11:31:23			
		EEG Trace	19.04.2021 10:46:53		EEG Video	23.04.2021 12:16:24			
		EEG Trace	19.04.2021 18:46:50		EEG Video	23.04.2021 13:01:24			
		EEG Trace	20.04.2021 2:46:46		EEG Video	23.04.2021 13:46:25			
		EEG Trace	20.04.2021 10:46:44		EEG Video	23.04.2021 14:31:25			
		EEG Trace	20.04.2021 18:46:41		EEG Video	23.04.2021 15:16:26			
		EEG Trace	21.04.2021 2:46:39		EEG Video	23.04.2021 16:01:27			
		EEG Trace	21.04.2021 10:46:38		EEG Video	23.04.2021 16:46:27			
		EEG Trace	21.04.2021 18:46:36		EEG Video	23.04.2021 17:31:28			
		EEG Trace	22.04.2021 2:46:33						
		EEG Trace	22.04.2021 10:46:30						
		EEG Trace	22.04.2021 18:46:28						
		EEG Trace	23.04.2021 2:46:25						
		EEG Trace	23.04.2021 10:46:22						


SystemPlus EVOLUTION | **Resource** | **Available GB** | **Descr**

Resource	Available GB	Descr
Archive2	345,26/2 047,40	
EEG	345,89/2 794,50	
Local Data	345,26/2 047,40	
Secure	345,26/2 047,40	
SQLbase	345,26/2 047,40	
SQLbase2	345,26/2 047,40	
Working Area	345,26/2 047,40	
Archive	345,26/2 047,40	1

Read/Write | Local Zone




582

SystemPLUS EVOLUTION


File Modify Acquire Analyze Balance Tools Window Help

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