
Manufacturer:

Siemens Healthcare GmbH
Henkestr. 127, 91052 Erlangen
Germany

User Manual supplement for Russia
from October 2022

***The magnetic-resonance system MAGNETOM Prisma
with accessories***

Friederike Hertle

Head of Regulatory Affairs
Siemens Healthcare GmbH, Magnetic Resonance (MR)

1 The name of the Medical Device

The magnetic - resonance system MAGNETOM Prisma with accessories

2 Manufacturing information of the Medical Device

2.1 Manufacturer

Name:	Siemens Healthcare GmbH
Address:	Henkestr. 127, 91052 Erlangen, Germany
Phone:	+49 9131 84-0
Web:	www.siemens-healthineers.com

2.2 Manufacturing site

Name:	Siemens Healthcare GmbH, Magnetic Resonance (MR)
Address:	Henkestr. 127, 91052 Erlangen, Germany
Phone:	+49 9131 84-0
Web:	www.siemens-healthineers.com

2.3 Local representative in Russia

Name:	Siemens Healthcare LLC
Address:	Dubininskaja Str. 96, 115093 Moscow, Russia
Tel.:	+7 495 737 12 52
Web:	www.siemens-healthineers.com/ru/

3 Designation of the medical device

3.1 Intended use

The MAGNETOM Prisma system is indicated for use as a magnetic resonance diagnostic device (MRDD) that produces transverse, sagittal, coronal and oblique cross sectional images, spectroscopic images and/or spectra, and that displays the internal structure and/or function of the head, body, or extremities. Other physical parameters derived from the images and/or spectra may also be produced. Depending on the region of interest, contrast agents may be used. These images and/or spectra and the

physical parameters derived from the images and/or spectra when interpreted by a trained physician yield information that may assist in diagnosis.

The MAGNETOM system may also be used for imaging during interventional procedures when performed with MR compatible devices such as in-room displays and MR Safe biopsy needles.

The MAGNETOM system is not a device with measuring function as defined in the EU directive 93/42/EEG (MDD) / EU regulation 2017/745 (MDR). Quantitative measured values obtained are for informational purposes and cannot be used as the only basis for diagnosis.

Your MR system is a medical device for human use only!

3.2 Indications

Magnetic Resonance Imaging (MRI) is a noninvasive technique used for diagnostic imaging. MRI with its soft tissue contrast capability enables the trained user to differentiate between various soft tissues, for example, fat, water, and muscle, but can also discriminate between the different bone structures. This allows for visualization of anatomic details and can offer more sensitivity and specificity in the detection of certain abnormalities than other imaging modalities, when read and interpreted by trained personnel. In general, it is the health care provider who decides on the appropriate diagnostic procedure based on, for example, the patient's physical examination and medical history.

3.3 Contraindications and possible side effects

This chapter contains very important information about MR safety: In principle, a qualified physician must evaluate the risk/benefit ratio of the MR examination for every patient.

To date, there is no scientific proof that MR examinations are harmless for pregnant women, the unborn (embryos or fetuses) and children under two years of age.

In general, an MR examination is contraindicated for patients with electronic or electronically conductive implants or metals, especially those containing ferromagnetic foreign matter.

Typical contraindications for MR examinations are:

- Electronic implants: e. g. pacemakers, stimulators, insulin pumps
- Artificial heart valves, aneurysm clips
- Metal splinters in the eye (danger of retinal detachment)
- Artificial anus (anus praeter) with magnetic closure
- Transdermal drug patches with metallic backings
- Electrically conductive implants and prostheses
- Metallic spirals for contraception (IUDs = Intrauterine Devices)

- Transdermal or other similar implants (for example, body piercings as well as magnetic piercings)

 WARNING
Electronic and/or electrically conductive implants and magnetizable inclusions in static and low-frequency magnetic fields and RF fields! Risk of patient death/patient injury ♦ Ask the patient about implants and inclusions. ♦ Do not perform MR examinations on patients with electronic or electrically conductive implants and magnetizable inclusions. ♦ Ensure that patients wearing such implants and/or inclusions remain outside the exclusion zone (0.5 mT line).

Exceptions: Certain implantable medical devices have been cleared, approved and/or licensed by the Competent Governmental Authorities and/or labeled by the device manufacturer as “MR-conditional”. For such implantable medical devices, the previously mentioned list of general contraindications and the warning may not be applicable in its entirety.

It is the responsibility of the device manufacturer to declare an implantable medical device as MR-conditional if appropriate and to define the conditions (constraints) for safe MR scanning. The MR operator must be aware of any such conditions for MR scanning. It is the obligation of the MR operator to assure that these conditions are strictly adhered to. To obtain these specific conditions the MR operator may refer to the labeling of the implantable medical device or contact the device manufacturer. Siemens MR does not assume responsibility or liability for the operation of the MR system with any implantable medical device. Especially, Siemens MR is not responsible for controlling technical parameters of the MR system other than those defined by the normal operating mode, the first level controlled operating mode and the data provided in the system owner manual, such as spatial gradient field plots.

 CAUTION
Eddy currents induced by low-frequency magnetic fields! Patient burns ♦ Do not examine patients with electrically conducting implants or prostheses.

CAUTION

Electrically conducting objects!

Injury to patient due to warming

Incorrect diagnosis due to artifacts

- ◆ Request that the patient removes all electrically conducting objects, e.g. necklaces, rings, braces, rubber bands for long hair, piercings as well as jewelry.
- ◆ Request that the patient removes all clothing including electrically conducting material, for example, bras, metallic appliques or woven metallic yarns.
- ◆ Inform patients that eyeliners and tattoos may contain ingredients causing artifacts or skin irritations during MR examinations. In some cases, patients have been burned.
- ◆ To prevent injuries, instruct patients to remove makeup prior to the examination.
- ◆ Instruct patients to seek medical attention in case of discomfort during or following the MR examination.

Possible side effects

Possible undesirable side effects for MR are dizziness, heating, claustrophobia and nerve stimulation.

4 Classification of the Medical Device

4.1 Type of the medical device in accordance with the nomenclature for classification

Europe (CE)

Kind of product:	Active medical device for diagnosis
Classification rules according to MDD 93/42/EEC, Annex IX:	Rule 10
Classification according to Annex IX:	Ila
Classification criteria:	Active devices intended for diagnosis are in Class Ila if they are intended to supply energy which will be absorbed by the human body, except for devices used to illuminate the patient's body, in the visible spectrum

Russia

Classification in accordance with Russia Nomenclature classification:	135160
Classification in accordance with Russian GOST R 50444-2020	Group 1

4.2 Class of the potential risk of the medical device application.

SW safety classification acc. to EN 62304

Based on the current product risk management results, the SW safety class of XA30 is "B".
The SW safety classification is done during risk management.

Classification of the Medical device according to MDD: IIa.

4.3 Electrical safety class:

System	Class I, in accordance with IEC 60601-1.
Gantry	Class I, in accordance with IEC 60601-1.
Patient table, in variants: Fixed; Dockable Tim table	Class I, in accordance with IEC 60601-1.
Physiological Measurement Unit (PMU)	Class I, in accordance with IEC 60601-1.

4.4 IP protection class

System	IP X0
Gantry	IP X0
Patient table, in variants: Fixed; Dockable Tim table	IP X0
Physiological Measurement Unit (PMU)	IP X0

4.5 Applied parts

Patient table	Type B applied part
Head/Neck 20 coil	Type BF applied part
Spine 32 coil	Type BF applied part
Body 18	Type BF applied part
Body 18 long	Type BF applied part
Flex Small/Large 4	Type BF applied part
UltraFlex small/large 18	Type BF applied part
Head/Neck 64	Type BF applied part
Tx/Rx CP Head Coil	Type BF applied part
Body 30/60	Type BF applied part
Peripheral Angio 36	Type BF applied part
TxRx 15 Ch Knee	Type BF applied part
TxRx Knee 15	Type BF applied part
TxRx Knee 15 Flare	Type BF applied part
Tx/Rx CP Extremity Coil	Type BF applied part
Hand/Wrist 16	Type BF applied part
Foot/Ankle 16	Type BF applied part
Head 32	Type BF applied part

Pediatric 16	Type BF applied part
Multi-nuclear spectroscopy coil Head Coil 13C1H Head Coil 31P1H Head Coil 23Na1H	Type B applied part
Shoulder 16 Coil Kit: Small/Large 16	Type BF applied part
8-ch. Sentinelle Breast Coil 16-ch Sentinelle Breast Coil	Type BF applied part
2-/4-/8-Channel Sentinelle Breast Coil 2-/10-/16-Channel Sentinelle Breast Coil.	Type BF applied part
Breast 18	Type BF applied part
4-channel BI Breast Coil 16Ch AI Breast coil	Type BF applied part
Endo Interface	Type BF applied part
Loop Coil, large	Type BF applied part
Loop Coil, medium	Type BF applied part
Loop Coil, small	Type BF applied part
Special Purpose Coil	Type BF applied part
Flex Coil Interface	Type BF applied part
Tim Coil Interface	Type BF applied part
Tim Coil Interface MNO	Type BF applied part
Multi-channel Interface	Type BF applied part
Body coil (integrated)	Type B applied part
ECG/Pulse module	Type BF applied part

4.6 Explosion protection

The MR system is not intended for operation in areas prone to explosion (e.g., highly flammable mixtures of anesthesia gases with air or oxygen or nitrous oxide)

4.7 Degree of safety in the presence of a flammable anesthetic's mixture with air or with oxygen or with nitrous oxide

No AP or APG category equipment

4.8 Operating mode

Continuous operation in accordance with p. 6.6 IEC 60601-1

4.9 Mains operated equipment with additional power sources

None

4.10 Laser light localizer

Laser class: 1M

Laser wavelength: 640 nm

Maximum output of laser radiation (accessible emission level, AEL): 1,41 mW

Used laser light localizer are Class 2M for older systems and Class 1M for newer systems. For the purpose of maintenance and service, the laser symbol label as well as the explanatory text label have to be permanently affixed to a location easily visible to the service personnel. The explanatory text must be provided either in local language or accepted alternative language depending on the country of installation.

The laser of the laser-light localizer is classified as Class 1M according to IEC 60825-1:2007 and IEC 60825-1:2014. Therefore, warning labels are not required.

5 Technical Data: Specifications and Characteristics

5.1 Design data

MAGNETOM Prisma with software *syngo* MR XA30 is a medical MRI system. The system contains advanced magnet technology including superconducting magnet technology. The large bore, whole body scanner is designed for increased patient comfort.

5.2 Photos of the device

Front view:


5.3 The list of components and accessories

No	Name of component	Description	Photos
	Basic block:		
1.	Whole-body superconductive 3T magnet with active shielding (AS) technology with counter coils	Whole-body superconductive 3T magnet with active shielding (AS) technology with counter coils.	
2.	Tim [204x64] XR Gradients #P Tim [204x128] XR Gradients #P	The whole-body gradient system XR engine	 

3.	Integrated no tune transmit/receive Body Coil	Integrated no tune transmit/receive Body Coil	
4.	Magnet cover	Magnet cover	
5.	Cooling system	Integrated magnet cooling system	

6.	Tim Table #P Tim Dockable Table #P (Add. Tim Dockable Table #P)	The new Tim Table is designed for maximized patient comfort and smooth patient preparation. The unique design of the Tim Table can support up to 250 kg (550 lbs) patients without restricting the vertical or horizontal movement. The Tim Dockable Table is designed for maximum patient comfort and smooth patient preparation. Tim Dockable Table can support up to 250 kg (550 lbs) patients without restricting the vertical or horizontal movement. The one step docking mechanism and the innovative multi-directional navigation wheel ensure easy maneuvering and handling. Critically ill or immobile patients can now be prepared outside the examination room for maximum patient care, flexibility and speed.	 
----	---	---	--


7.	Electronics cabinets: - Gradient cabinet - Control cabinet	Electronics cabinets: - Gradient cabinet - Control cabinet	
----	--	--	--

8.	Separator 60kW	Function: Interface between the on-site water chiller (of any brand/type) or Interface to the central hospital chilled water supply.	
9.	The alarm box	The alarm box has the following functions: - Displays alarm signals - Switches the MR system on and off - Magnet Stop/ Magnet Quench	 


			
--	--	--	--

			
10.	The intercom	The intercom allows personnel and patients to communicate during the examination	


11.	syngo Acquisition Workplace	syngo Acquisition Workplace is the workplace in the control room. It includes the host processor with the operating elements monitor and mouse.	 
12.	PC Keyboard US English #NX	The MR system is equipped with an original Siemens keyboard. This keyboard is a modified Windows keyboard where the numeric keys have been replaced with symbol keys.	
13.	Head/Neck 20 coil	Head/Neck 20 Coil provides examinations of the head and the neck (e.g. nape of the neck and vessels of the neck), and combined measurements with other coils	
14.	Spine 32 coil	The Spine Coil is used for examining the spine and for combined measurements with other coils	


15.	Body 18 Body 18 long	The Body Coil enables examinations of the thorax, abdomen, pelvis or hip. It can also be used for cardiac or vascular applications.	
16.	Flex Coil Interface 3T	Interface with integrated preamplifiers for the connection of the following flexible coils to the system: - Flex Large 4 - Flex Small 4 - Loop Coil, large - Loop Coil, medium - Loop Coil, small The interface is not permanently mounted and therefore allows free positioning of the flexible coils as required by the examination procedure.	
17.	Flex Small/Large 4 or UltraFlex small/large 18	The Flex coil is available in sizes small and large. Depending on the area to be examined, the coil is positioned on the patient and is used to image the area of the body it covers.	Flex Small/Large 4 : 

			 UltraFlex small/large 18:  
--	--	--	---


18.	Head/Neck 64	The Head/Neck 64 can be used for applications like head examinations, neck examinations, MR Angiography, combined head/neck examinations or for imaging of the TMJ (temporo-mandibular joints).	
19.	Tx/Rx CP Head Coil	Tx/Rx CP Head Coil is suitable for the following examinations: - Vascular examinations in the head - E. g., vascular examinations of the base of the skull and the cerebrum	
20.	Body 30/60	The coil with its integrated pre-amplifiers ensures excellent signal-to-noise ratio. The light-weight coil ensures highest patient comfort. The highly flexible design enables a wide variety of applications including: - Thorax (incl. heart) - Abdomen - Pelvis - Hip The Body 30 coil allows you to examine the following parts of the body:	Body 30/60 


		■ Upper part, stomach, liver ■ Pelvis ■ Heart Depending on the requirements of the examination, the coil may be used in the Head First or Feet First orientation. The Body 30 coil can be used as a front coil in conjunction with a spinal coil, or in conjunction with a second Body 30 coil as a rear coil. The combination of two Body 30 coils (front/rear) is called "Body 60".	
21.	Body 18 -> Body 18 long	This option exchanges the Body 18 from the standard coil configuration for the Body 18 long. The long cable of the Body 18 enables flexible combinations with additional Body 18, the Body 30 or Body 60 for extended anterior coverage.	

22.	Peripheral Angio 36	The Peripheral Angio 36 coil is suitable for examining lower extremity vessels. Depending on the requirements of the examination, the coil may be used in the Head First or Feet First orientation. The Peripheral Angio 36 coil can be used with a spinal coil and one or more body coils. By combining these coils, the examination area expands.	
23.	TxRx 15 Ch Knee TxRx Knee 15 TxRx Knee 15 Flare	The TxRx 15 Ch Knee coil is used for examinations of the left or right knee. The TxRx Knee 15 coil is also used for examinations of knees. The coil can be moved on the base plate along the longitudinal axis to bring it to the center to obtain an image of optimal quality. The TxRx Knee 15 Flare Knee Coil is used to examine the left or right knee joint. The coil can be moved on the base plate along the longitudinal axis to bring it to the center to obtain an image of optimal quality. New 15-channel transmitter/receiver coil for joint examinations in the area of the lower extremities.	TxRx 15 Ch Knee  TxRx Knee 15 TxRx Knee 15 Flare 



		Main features : - 15-element design (3x5 coil elements) with 15 integrated preamplifiers - iPAT-compatible - SlideConnect Technology	
24.	Tx/Rx CP Extremity Coil	Circularly Polarized no-tune transmit/receive coil for joint examinations in the region of the lower extremities.	
25.	Hand/Wrist 16	The Hand/Wrist 16 coil is intended for examinations of wrists, hands, fingers. Depending on the requirements for the examination, the coil can be used on a patient with the Head First orientation (position on the stomach) or Feet First orientation (position on the back). For Feet First orientation, the coil is mounted on the base plate.	

26.	Foot/Ankle 16	Foot/Ankle 16 is used for examining the feet and ankles.	
27.	Head 32	Head 32 provides: - High resolution head proton imaging - MR angiography of the head - Functional imaging of the brain	
28.	Pediatric 16	New 16-channel receive coil for head and neck imaging of newborns and children up to 18 months of age. Main Features: - 16 elements, 13 elements in the head and 3 elements in the neck - iPAT-compatible - Can be used in combination with all other Tim coils	 

		- Open anterior design - Infant cradle	
29.	Multi-nuclear spectroscopy coil Head Coil $^{13}\text{C}^1\text{H}$ Head Coil $^{31}\text{P}^1\text{H}$ Head Coil $^{23}\text{Na}^1\text{H}$	Non-proton nuclei, such as ^{31}P, ^{23}Na and ^{13}C can provide important information about the metabolism of the brain and diseases of the brain including stroke, tumors and inflammation. Dual tuned volume resonator with open design: - ^{31}P transmit / receive - quadrature polarization - resonance frequency 49.9 MHz - ^1H transmit / receive - quadrature polarization - resonance frequency 123.2 MHz Dimensions: • Head Coil $^{23}\text{Na}^1\text{H}$ - measurements 22 x 33 cm • Head Coil $^{13}\text{C}^1\text{H}$ - length 24 cm; outer diameter 35 cm; inner diameter 26.5 cm • Head Coil $^{31}\text{P}^1\text{H}$ - length 24 cm; outer diameter 35 cm; inner diameter 26.5 cm	Head Coil $^{13}\text{C}^1\text{H}$ and Head Coil $^{31}\text{P}^1\text{H}$  Head Coil $^{23}\text{Na}^1\text{H}$ 
30.	Shoulder 16 Coil Kit: Small/Large 16	The iPAT compatible Shoulder 16 Large and Shoulder 16 Small are ergonomically designed and adapted to the shape of the shoulder.	

		The different sizes obtain maximum image quality for different body sizes: - 165 mm (6.5 in) diameter for small and medium sized shoulders - 200 mm (7.9 in) diameter for large shoulders The coils can be used either for left or right shoulders. It features sliding attachments to the base plate and can easily be adjusted for comfortable positioning. The coils excels in highest resolution imaging with exceptional signal/noise ratio.	
31.	Coils for mammary glands with a modular design - no more than 2 pcs. (if necessary), versions: - 8-ch. Sentinelle Breast Coil; - and/or - 16-ch Sentinelle Breast Coil.	The breast coil with modular design delivers brilliant image quality for high-resolution 2D and 3D MR breast imaging. Techniques for reducing scan times, such as parallel imaging, can be used as well.	8-ch. Sentinelle Breast  16-ch Sentinelle Breast Coil 


32.	Coils for mammary glands with a modular design for imaging and/or biopsy - not more than 2 pcs. (if necessary), versions: - 2-/4-/8-Channel Sentinelle Breast Coil; and/or - 2-/10-/16-Channel Sentinelle Breast Coil.	Coils for mammary glands with a modular design for imaging and/or biopsy is used for patient examinations, for biopsy, or for conventional imaging of the mammary gland.	
33.	Breast 18	The Breast 18 is used for examining female patients. Depending on the requirements of the examination, the coil may be used in the Head First or Feet First orientation. Both sides of the breast are displayed simultaneously to allow for comparative evaluation	 
34.	4-channel BI Breast Coil 16Ch AI Breast coil	The 4-channel BI Breast Coil is well suited for breast biopsy with lateral, medial and cranio-caudal access.	4-channel BI Breast Coil

			 16Ch Breast Coil 
35.	Endo Interface	Kit including interface device for connecting the prostate receive coil probe. The interface device connects to the Tim Coil Interface.	
36.	Tim Coil Interface	Coil adapter plug for up to 8 receive and 1 transmit channels, in order to connect existing dedicated knee and breast coils (Tx/Rx 15-channel Knee Coil, 4-channel BI Breast Coil, 16-channel AI Breast Coil, (2/4)/8-channel Sentinelle BreastCoil and (2/10)/16-channel Sentinelle BreastCoil)	

37.	Tim Coil Interface MNO	Coil adapter plug for up to 8 receive and 1 transmit channels, in order to connect dedicated multinuclear coils with MAGNETOM Prisma.	
38.	Multi-channel Interface #3T	The Multi-Channel Interface transforms the Head Coil Direct Connect sockets into two Tim 4G Slide Connect sockets. In combination with Tim Coil Interfaces (not included in this option), this allows to connect specific Legacy or 3rd party coils.	
39.	UltraFlex Large 18 Flex Large 4	Depending on the area being examined, the coil is mounted on the patient and allows you to display the area of the body that it closes. For example, you can examine the following parts of the body: ■ Hand ■ Elbow and shoulder ■ Hip and knee	Flex Large 4  UltraFlex Large 18

			
40.	UltraFlex Small 18 Flex Small 4	Depending on the area being examined, the coil is mounted on the patient and allows you to display the area of the body that it closes. For example, you can examine the following parts of the body: ■ Hand ■ Elbow and shoulder ■ Hip and knee	UltraFlex Small 18  Flex Small 4 


41.	Loop Coil, large	No-tune loop receive coil for examinations of the upper or lower extremities (e.g. shoulder).	
42.	Loop Coil, medium	No-tune loop receive coil for examinations of small structures near the surface (e.g. hands, fingers, skin).	
43.	Loop Coil, small	No-tune loop receive coil for examinations of small structures near the surface (e.g. hands, fingers, skin).	
44.	Special Purpose Coil	The special purpose coil is a no-tune receiver coil designed for examination with small Field-of-Views.	
45.	Connection cables	Connection cables	




			
46.	Phantoms	Measurement phantoms are used for quality measurements.	Plastic bottle 15ml:  Plastic bottle 1900ml:  Plastic bottle 5300ml:  Spherical phantom D165:  Spherical phantom D240: 

			 Hand/Wrist 16 phantom:  Foot/Ankle 16 phantom:  Special phantom for Endorectal coil: 
47.	Phantom holders	Phantom holder	Phantom holder 17:  Phantom holder 18

			 Phantom holder 31
			 Phantom holder 32
			 Phantom holder 25
			 Phantom holder 26
			

			Phantom holder 27  Phantom holder 28  Holder for spheric phantom  Phantom holder 29  Phantom holder 40  Phantom holder 45
--	--	--	---


			 Phantom holder Breast Sen  Endorectal coil simulator 
48.	The application package "Tim Application Suite": -Basic software license: basic license key for running the syngo software; -Basic MR applications license; -3D functionality; -Argus (cardiac) basic viewing; -DICOM worklist support; -PhoenixZIP; -Different standard measurement/examination and post-processing functionality.	The Tim Application Suite offers a complete range of clinically optimized sequences, protocols and workflow functionalities for all body regions.	

49.	ECG electrodes, non-magnetic, disposable (if necessary, no more than 2000 pieces)	Disposable electrodes for long-term monitoring For children and adults X-ray transparent	
50.	Additional PMU Sensor Kit #T+D Additional PMU Sensor Kit #T+D	Physiological Measurement Unit (PMU) - a device for monitoring physiological parameters with wireless control and triggering; synchronizes the measurement with the physiological cycles of cardiac and/or respiratory movements.	 ECG/Respiratory Unit  Peripheral Pulse Unit 

51.	MRI-Stretcher, non-magnetic;	MRI-Stretcher, non magnetic, chrome-nickel compounded for patient transport inside the clinic building. Continuously adjustable head piece with spring brake, synthetic leather cushion, movable barrier on either side, stretcher movable on four fixable wheels diameter 100 mm (for max. 3 mm floor unevenness). Dimensions: 190 x 65 x 80 cm Maximum laminary load: 200 kg Weight (empty): 41 kg Heights of 65, 70, 75 and 80 cm are available without additional costs. Please send an info with the order which height is favoured. Standard upholstery colour: blue	
-----	------------------------------	--	---


52.	Operator Manual – System and Data Management	Instruction for Use	 The cover features a large red header. Below it, the text 'syngo MR XA30' is displayed in red, followed by 'Руководство пользователя – Управление системой и данными' in black. At the bottom right, the Siemens Healthineers logo is visible.
53.	Operator Manual – TimTX TrueShape	Instruction for Use	 The cover features a large red header. Below it, the text 'syngo MR XA30' is displayed in red, followed by 'Руководство пользователя – TimTX TrueShape' in black. At the bottom right, the Siemens Healthineers logo is visible.

54.	Operator Manual – Vascular MR	Instruction for Use	
55.	Operator Manual – Orthopedic MR	Instruction for Use	

56.	Operator Manual – Scanning	Instruction for Use	
57.	Operator Manual – Coils	Instruction for Use	


58.	Operator Manual – Dot Cockpit	Instruction for Use	
59.	Operator Manual – IT Configuration and Security	Instruction for Use	

60.	Operator Manual – MR System	Instruction for Use	
61.	Operator Manual – Neuro MR	Instruction for Use	

62.	Operator Manual – Oncologic MR	Instruction for Use	
63.	System Owner Manual	Instruction for Use	

64.	Quick Guide	Instruction for Use	
65.	Operator Manual – Body MR	Instruction for Use	

66.	Operator Manual – Breast MR	Instruction for Use	
67.	Operator Manual – Cardiac MR	Instruction for Use	

68.	Operator Manual – 2-/4-/8-Channel Sentinelle Breast Coil, -/10-/16-Channel Sentinelle Breast Coil	Instruction for Use	
69.	Operator Manual – 4 Ch BI Breast / 16 Ch AI Breast	Instruction for Use	

70.	Addendum – Feature List MAGNETOM Family	Instruction for Use	
71.	Addendum – PMU information MAGNETOM Family	Instruction for Use	

72.	Addendum – Regulatory Information MAGNETOM Family	Instruction for Use	
73.	Addendum – Use of syngo.MR in Numaris/X	Instruction for Use	

74.	Operator Manual – Access-i	Instruction for Use	
75.	Operator Manual – Advanced User	Instruction for Use	


76.	Operator Manual – Measurement Parameters	Instruction for Use	
77.	Operator Manual – MR Fingerprinting	Instruction for Use	


78.	Operator Manual – MR View&GO	Instruction for Use	
79.	Release Information	Instruction for Use	


80.	User manual supplement for Russia	Instruction for Use	
Accessories:			
1.	syngo MR Workplace #NX PC Keyboard US English #NX	The syngo MR Workplace allows for evaluation, documentation and post-processing of previously measured images while acquiring images at the syngo Acquisition Workplace.	

			 
2.	Workstation-based One User Solution	syngo.via Workstation is a syngo.via-based entry level sales package for advanced visualization customers. It provides the same software functionality, as the conventional syngo.via product.	
3.	PC Keyboard US English	The MR system is equipped with an original Siemens keyboard. This keyboard is a modified Windows keyboard where the numeric keys have been replaced with symbol keys.	

4.	Foot switch In-Room #T+D	MR Foot Switch In-Room with two pneumatic buttons for start and stop of a preset MR sequence. The foot switch is MR compatible and is positioned near the patient table on the examination room floor. The MR foot switch is especially suited for sequence control when MR examinations are carried out during operations.	
5.	In-Room syngo Acquisition Workplace	The In-Room syngo Acquisition Workplace is an additional operating console in the examination room. It is used for image viewing and MR system operation and comprises the following components: - LCD monitor - Tray with trackball and keys - Foot switch (optional)	

6.	Supervision	The supervision monitors that cables, connections, or other components of the gradient system do not show excessive heating.	
7.	Respiratory cushion with pressure hose	Respiratory cushion with pressure hose	
8.	Respiratory belt	Respiratory belt	
9.	Application cushion	Cushion Peru	

10.	Charging station for PMU	The charging station is installed separately near the syngo Acquisition Workplace and is used for storing the PERU (Physiologic ECG and Respiratory Unit) and PPU (peripheral pulse unit) that are part of the Physiological Measurement Unit.	
11.	Add. Dot Control Center Rear #P Dot Display	The Dot (Day Optimized Throughput) display provides status information as well as several tools. The display is located above the magnet bore opening at the front of the magnet cover. Optionally, at the back of the magnet, a second display can be installed for the additional control unit.	
12.	Squeeze bulb	Patients may use the squeeze bulb to alert the operating personnel.	


13.	Headphones	The headphones are intended for communication with the patient and can be used in combination with ear plugs.	
14.	MagnaCoil Headset	The MagnaCoil Headset can be used as an alternative headphone (for example, in head coils) to the standard Siemens headphones which are part of the products since years. Sometimes the standard headphones are not very practicable and therefore Siemens wants to offer an alternative: the MagnaCoil Headset.	
15.	Adapter plate for Tx/Rx CP Extremity coil +TxRx Knee coil	Adapter plate for connecting the coils to MAGNETOM Prisma	


16.	Accessory Breast 18 60cm Accessory Breast 18 70cm	The accessories kit for the Breast Coil contains a table pad and the connection wedge. The connection wedge is used as an abdominal support cushion for the patient.	
17.	Positioning frame f. BreastCoil	Open design positioning frame for inserting the 4-channel BI Breast Coil or 16-Channel AI Breast Coil. The positioning frame comes supplied with: - a head cushion - fixation rail to proper position and fix the coil positioning frame on the patient table - abdominal support cushions - Breast stabilizing cushions - positioning guidance	Frame: (white part of the coil): 
18.	Biopsy Starter Kit #Tim	The breast biopsy kit comes as a case containing a fixation unit plus accessories for both the Post-Pillar and the Grid biopsy method. Please note that the scope of delivery from Siemens provides biopsy accessories only for training purposes. Accessories for the examination (needles, needle guide, height adjustable grid, etc.) need to be ordered separately, for example, from NORAS MRI products GmbH.	


19.	Cushion Set 60cm #T+D	The Cushion Set 60cm is required for adapting Sentinelle Breast Coils for use with 60 cm Tim 4G systems.	
20.	Vacuum cushions, small	Vacuum cushions for stable positioning of the patient during the examination. The kit comprises of an integrated vacuum pump and 3 anatomically shaped (spine, head, small joints) cushions of different sizes. The Comfort Kit supports a more comfortable positioning of the patient and improves patient experience by helping to reduce pain.	Small: 
21.	Vacuum cushions head	Vacuum cushions for stable positioning of the patient during the examination. The kit comprises of an integrated vacuum pump and 3 anatomically shaped (spine, head, small joints) cushions of different sizes. The Comfort Kit supports a more comfortable positioning of the patient and improves patient experience by helping to reduce pain.	Head: 
22.	Vacuum cushions spine	Vacuum cushions for stable positioning of the patient during the examination. The kit comprises of an integrated vacuum pump	Spine:

		and 3 anatomically shaped (spine, head, small joints) cushions of different sizes. The Comfort Kit supports a more comfortable positioning of the patient and improves patient experience by helping to reduce pain.	
23.	Trigger Converter	An optical trigger signal is available to trigger external stimulation devices in fMRI experiments.	
24.	Patient Supervision TV Supervision Camera &Controller #T+D Additional Camera #T+D	This package contains a special video camera for monitoring the patient during an MR examination. The supervision solution is customizable and designed to address different site specific requirements. Up to 4 cameras can be optionally connected for patient supervision in the examination or waiting room.	

		Additional video camera to allow for comprehensive monitoring of the patient's face independent of patient position during an MR examination.	
25.	Patient TV wall support Patient TV wall support (multi-cam)	Wall mount for the patient monitor.	
26.	LCD Color Remote Monitor #NX	This remote monitor is used for connection in parallel to the syngo Acquisition Workplace monitor. The LCD remote monitor is of the same type as the LCD color monitor of the syngo Acquisition Workplace with the same high-quality characteristics: - can be tilted horizontally forward and backward, - integrated gamma correction for optimum display of radiographic grayscale, automatic backlight control for long-term brightness stability.	
27.	Video Distributor for monitor #T+D	Video distributor for connection of up to two additional color LCD monitors as well as an In-Room Acq WP syngo (in the examination room parallel to the syngo Acquisition Workplace operator monitor).	


28.	Coil Storage Cart #T+D,Ez	Specially designed non-ferromagnetic cart for easy storage of the most commonly used coils and accessories. The cart may be rolled to convenient locations in the examination room and can be opened up to work like a shelf. The coil storage cart has multiple drawers and trays as well as many other storage spaces for coils, cushions and miscellaneous items. Its dimensions are: Width 140 cm (4' 7") when closed and 280 cm (9' 12") when opened, depth 54 cm (1'9") and height 121 cm (3'12").	 Coil storage cart Peripheral Angio 36 
29.	MR coil shelf	HAMACO accessorial shelf inside the rf-shielded enclosure (coil shelf). L=2.00 x W=2.10 x D=0.50 m Weight: 250 kg Same design as the rf-enclosure inner finish with 2 plates 1 x 1.20, 1 x 0.80 m and 5 pieces enforced plates adjustable in the vertical position, with additional 45° end plates. The shelf is pre-assembled and fragmented as kit available.	



30.	Accessories cabinet	Specially designed non-ferromagnetic cabinet for easy storage of some of the most commonly used coils and accessories. - Width 913 mm - Depth 725 mm - Height 2050 mm Weight 214 kg	
31.	Cushions and positioning aids - table pad 1 - table pad 2 - table pad 3 - table pad 4 - head-heel insert - 4 belt straps - different-sized sandbags - lookout mirror holder	Accessories for patient positioning during examination	Table pads, available in four different sizes:  Head-heel insert: 

			 Belts and straps:   Different-sized sandbags:  Mirror holder with mirror:
--	--	--	--

			
32.	Set of positioning aids - knee wedge - positioning roll - back-of-the-head support - large triangles - small triangles - wedge - lordosis cushion - flat cushions - flat wedges - knee cushions - different-sized multipurpose cushions - RF homogeneity cushion	Accessories for patient positioning during examination	Positioning roll:  back-of-the-head support:  Knee support:  Knee cushion:


			 Triangles small/large:  Wedge:  Flat wedges:  Lordosis cushion:
--	--	--	---



			 Multipurpose cushions:  Flat cushions:  RF homogeneity cushion: 
33.	Other accessories: - paper roll holder - protective film for tabletop	Designed for patient positioning.	Paperroll holder:  protective film for tabletop:


			
34.	Head/Neck 20 accessories: - pad 32 - Head/neck cushion - head stabilization wedges - clamp cushions - mirror carrier with mirror - protective cover for the head/neck coil sockets	Designed for patient positioning.	Pad 32:  Head/Neck cushion:  Head stabilization wedges:  Clamp cushions:


			 Mirror carrier with mirror:  protective cover for the head/neck coil sockets:  Head/neck 64 relax cushion  Head positioning wedges
35.	Head/Neck 64 accessories: - Head/neck 64 relax cushion- head positioning wedges - mirror carrier with mirror - clamp cushions	Designed for patient positioning.	



			 mirror carrier with mirror  clamp cushions 
36.	Head 32 accessories - pad 32 - Head/neck cushion - head stabilization wedges - mirror carrier with mirror	Designed for patient positioning.	Pad 32  Head/Neck cushion:






			 Head stabilization wedges:  Mirror carrier with mirror: 
37.	TxRx CP Head accessories - Head cushion - head wedges (thin size) - Head-Spine table pad - head wedges (thick size)	Designed for patient positioning.	Head cushion:  Head wedge thin:



			 Head-Spine table pad:  Head wedge thick: 
38.	Pediatric 16 accessories: (2) Pediatric cushions - head/crown pad - neck elevation pad - coil transition pad (for the use without infant cradle) (4) infant cradle (3) 2 straps (for infant cradle only) - cradle pad (for infant cradle only)	Designed for patient positioning.	 Pediatric 16 accessories, where 1 - head/crown pad 2 - neck elevation pad 3 - infant cradle 4 - 2 straps (for infant cradle only)


			5 - cradle pad (for infant cradle only)  Head/crown pad  Neck elevation pad  coil transition pad (for the use without infant cradle)  cradle pad (for infant cradle only)
--	--	--	--


			 infant cradle  2 straps (for infant cradle only) 
39.	Spine accessories: - Spine support	Designed for patient positioning.	Spine support 
40.	Body 18 accessories - straps - belts - Cardiac spacer cushion	Designed for patient positioning.	Belt:



			 Strap:  Cardiac spacer: 
41.	Body 30/60 accessories: - table pad 1 for Body 30 - table pad 2 for Body 30 - table pad 3 for Body 30 - table pad 4 for Body 30 - 2 belts 1500 - 4 belts Body 30 - 4 belt clips - frame for use as a rear coil;	Designed for patient positioning.	



	- socket adapter for Spine coil.		     
--	---	--	--



			
42.	PA 36 accessories -PA support	Designed for patient positioning.	
43.	Foot/Ankle 16 accessories - Knee cushion - Heel cushion - Foot cushion	Designed for patient positioning.	Knee cushion:  Heel cushion: 


			Foot cushion: 
44.	Hand/Wrist 16 accessories - Wrist pad - Elbow pad - Palm pad - Wedge cushion - Wrist support cushion	Designed for patient positioning.	Wrist pad:  Elbow pad:  Palm pad:  Wedge cushion:  Wrist support cushion:



			
45.	Flex Small/Large 4/Loop accessories - Hand wedge - soft straps - belts - iPat cushion, large - iPat cushion, small - cable holders (3 pcs.)- Ultraflex positioning pad - UltraFlex foot/ankle pad - soft wedge pad	Designed for patient positioning.	Hand wedge:  Soft straps:  Belts:  iPAT cushion small:



			 iPAT cushion large:  cable holders  Ultraflex positioning pad  UltraFlex foot/ankle pad 
--	--	--	---



			soft wedge pad 
46.	Endorectal accessories - Endorectal Adapter 3T - Endo cushion	Designed for patient positioning.	Endorectal Adapter 3T:  Endo cushion: 
47.	Breast 18 accessories: - coil pad - slider pads - body cushion 60 (wedge connector 60) - body cushion 70 (wedge connector 70) - head support - head support cushion - armrest pad - cushion for the breast fixations - table pad Breast 18	Designed for patient positioning.	 

			     
--	--	--	--





48.	Shoulder Shape 16 accessories: - L-shaped holder	Designed for patient positioning.	
49.	Shoulder Small/Large 16 accessories: - body pad - shoulder cup cushion	Designed for patient positioning.	
50.	TxRx 15 Ch Knee/TxRx CP Extremity accessories - Coil pad, 0.25 - Coil pad, 0.50 - Ankle pad - Knee rest pad - Foot rest pad	Designed for patient positioning.	Coil pad different sizes:  Ankle pad:  Knee rest pad: 


			 Foot rest pad: 
51.	TxRx Knee 15 Flare accessories - Coil pad, 0.25 - Coil pad, 0.50 - Ankle pad - Knee rest pad - Foot pad	Designed for patient positioning.	Coil pad different sizes:  Ankle pad:  Knee rest pad: 


			Foot pad: 
52.	TxRx 15 Ch Knee TxRx CP Extremity accessories Coil pad, 0.25 Coil pad, 0.50 Knee cushion Foot ankle rest pad Foot support 098	Designed for patient positioning.	    


53.	4Ch Breast BI, 16Ch AI Breast coils accessories - head rest; - abdominal pad; - foot support roller.	Designed for patient positioning.	
54.	16Ch AI Breast coils accessories Breast coil (for expert studies of mammary glands) (16 Ch AI Breast) accessories (no more than 30 pcs., if necessary): - Distance cushion; - Cushion to stabilize the breast; - Support mat; - Medial cushion; - Wing cushion large; - Wing cushion small:	Designed for patient positioning.	




55.	SpectroShim Support #P	The SpectroShim option provides additional shimming power to further improve inhomogeneities caused by the patient, volunteer, or measurement object.	
56.	2/4/8-ch Sentinelle Breast accessories 2/10/16-ch Sentinelle BreastCoil accessories - medial plug; - head rest; - body tabletop pad; - body wedge pad; - head rest pad; - plastic Cover for 4 Ch coil element; - shoulder wedges; - foot rest pad; - front pad; - hip wedge pad; - insulation pad; - shoulder support.	Designed for patient positioning.	


57.	2/4/8-ch Sentinelle Breast biopsy accessories 2/10/16-ch Sentinelle BreastCoil biopsy accessories - fixed marker block/ Marker block; - contra lateral support (for medial access); - immobilization plate; - grid plate.	- Fixed marker block/Marker block - Contra lateral support (for medial access) - Immobilization plate - Grid plate Please note that the scope of delivery from Siemens provides biopsy accessories only for training purposes. Accessories for the examination (needles, needle guide, height adjustable grid, etc.) need to be ordered separately, for example, from NORAS MRI products GmbH.	   
58.	-Optional Dot Engine Workflows;	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses,	

		depended on the individual customer order. • DotGO Routine Package #NX • DotGO Advanced Package #NX • DotGO XL Package #NX • Spine Dot Engine #T+D • Large Joint Dot Engine • Abdomen Dot Engine • Prostate Dot Engine • Cardiac Dot Engine #T+D • Angio Dot Engine #NX • Breast Dot Engine #T+D • LiverLab Dot Engine XA30/XA31 • Whole-Body Dot Engine #NX	
59.	-Optional Application Packages;	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order. • Kinetic Sensor Motion Correction • FREEZEit+ • Compressed Sensing: GRASP-VIBE, SPACE, TOF, SEMAC • LiverLab • MapIt • SWI (Susceptibility Weighted Imaging) • Inline BOLD Imaging • Simultaneous Multi-Slice (SMS) • Inline Perfusion • 3D PACE • DTI (Diffusion Tensor	

		Imaging) • DSI (Diffusion Spectrum Imaging) • TWIST • NATIVE • QISS • ASL (Arterial Spin Labeling): 2D, 3D • Compressed Sensing Cardiac Cine • MyoMaps • Flow Quantification • Advanced Interactive Realtime • RT respiratory self-gating • Single Voxel Spectroscopy • Chemical Shift Imaging for CSI 2D, CSI 3D • Multinuclear Support • Multinuclear Spectroscopy • MR Elastography • Breast Biopsy Software • IDEA (Integrated Development Environment for Applications)	
60.	-Optional Postprocessing Packages.	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order. • syngo.MR Onco • syngo.MR 3D Lesion Segmentation • syngo.MR OncoTrend • syngo.MR Onco Engine • syngo.MR OncoCare • syngo.MR Tissue 4D • syngo.MR BreVis	

		• syngo.MR Neuro Perfusion • syngo.MR Neuro Perfusion Mismatch • syngo.MR Brain Morphometry • syngo.MR Neuro fMRI • syngo.MR Tractography • syngo.MR Neuro Perfusion Engine • syngo.MR Spectro SVS • syngo.MR Spectro CSI • syngo.MR Spectro Engine • syngo.MR Spectro Extension • syngo.MR Spectro Research • syngo.MR Vascular Analysis • syngo.MR Cardiac 4D Ventricular Function • syngo.MR Cardiac Flow • syngo.MR Cardiac Perfusion • syngo.MR Cardio Engine • syngo.MR Cinematic VRT • syngo.MR Anatomic Intelligence	
61.	-Additional optional single Applications;	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.	


62.	-Argus Flow -Argus Dynamic Sygnal	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.	
63.	syngo Expert-i XA30	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order.	
64.	Syngo Security package	Software, optional application package. Software option for general regulatory security rules, providing enhanced security features including user management and audit trail functionality. This package supports customers in their achieving compliance with the HIPAA "Privacy" rule.	
65.	MR Workplace Table 1.2m MR Workplace Table, height adjust.	Table suited for syngo Acquisition Workplace and	

	Control Room Desk 1.35m	syngo MR Workplace based on syngo Hardware.	
		Workplace for operating personnel made of wooden, plastics and aluminum parts. Includes mount for electronic components.	
66.	MR Workplace Container, 50cm	50 cm wide extra case for the syngo host computer with sliding front door to allow change of storage media (CD/DVD/USB).	
67.	Adapter, RF Cabin #Tim,Ez	2 adapter plates for adjusting the larger filter plate opening of an existing RF-room enclosure. The new RF filter plate of a MAGNETOM MR system (current version) is incorporated into the same opening.	Filter plate with cabling: 


68.	RF-Cabin, Steel #3T	Steel full-room RF cabin, module method of building, for MAGNETOM Trio a Tim System, Verio and Prisma. Standard cabin outside dimensions: Maximum outside dimensions: Length x Width up to 35 sqm Height up to 3,00m Screen damping value: 90 dB (Co-Siting 100 dB) at 15 until 130 MHz.	 
69.	RF-Traverse Beams, Steel #3T	5 traverse beams for the RF cabin (Art. No. 14413606 and 14441589) of the MAGNETOM Trio a Tim System, Verio, Prisma as well as Biograph mMR. Length 4,590 mm	Lifting device without beams:

			 Beams, to use with lifting device 
70.	RF-Interior, LED #3T	Standard interior décor with LED lighting for the RF cabin (Art. No. 14413606 and 14441589) of the MAGNETOM Trio a Tim System, Verio, Prisma as well as Biograph mMR. Maximum outside dimensions: Length x Width up to 35 sqm Height up to 3,00 m	

71.	Filter Cabinet, big	Cabinet for filter plate inside of the RF cabin. Outside dimensions appr. 1.300 x 420 x 2.700 mm (L x W x H)	
72.	RF Cabin Sound Proofing, Steel #3T	Sound protection system for MAGNETOM Trio a Tim System, Verio and Prisma steel modular cabin (Art. No. 14413606) as well as Biograph mMR steel modular cabin (Art. Nr. 14441589). Maximum outside dimensions: Length x Width up to 35 sqm Height up to 3,00 m Due to building requirements according to DIN 4109 (sound protection for structural engineering) it is possible to get the following values: Sound protection, walls > 52 dB(A) Sound protection, ceiling > 56 dB(A) Sound protection, observation window > 50 dB(A)	Sound protection RF-Cabin: 

73.	Seismic Kit #T+D	Earthquake Kit for ECO Chiller 60 kW.	
-----	------------------	---------------------------------------	--


			
74.	Eco Chiller 60kW	The KKT ECO 133 - L chiller is a dedicated 20°C cooling system which automatically adapts to the different cooling requirements (e.g. system in operation, standby, ...) to reduce the energy consumption for cooling. The cooling system must be used in combination with the IFP (Interface Panel), if there is no on-site chilled water supply at all. The IFP is included in the scope of supply.	  



			
75.	High-End Computing [204x64] #NX High-End Computing [204x128] #NX	syngo MR XA30 hardware upgrade option, containing a new high performance measurement and reconstruction system.	
76.	Dual Monitor Package XA30	The Dual Monitor Package provides a second 24" LCD monitor for the acquisition workplace, identical to the system main host monitor. The two monitors provide space for protocol planning and exam progress on the left monitor, as well as viewing and post-processing functionalities on the right monitor. The Dot Cockpit can be used on both monitors as a floating window. This improves the MR examination workflow	

		by a smoother and more comfortable work space that avoids interruptions between planning, scanning, viewing and post-processing. It allows to keep running patient examinations always in sight to allow for fast interactions.	
77.	Advanced Host Computer #NX	The Advanced Host Computer offers increased computing power and increased memory for supporting an external syngo MR Workplace (optional) and/or to give a performance boost to applications that generate and process large data sets (e.g. BOLD imaging, fMRI post-processing).	


78.	UPS system with UPS Battery module	UPS system Liebert GXT4 3000RT230E for safeguarding computers. Including Power Cable of 9 m for connecting the UPS. UPS battery module Liebert GXT4 72VBATTE for safeguarding computers.	
79.	Filter plate for RF-Cabin	Filter plate in low leakage circuitry for shielded rooms (MRI) 6 separate lines, 250 V, 50/60 Hz, 250Vdc, 32A For usage at MR-Shielding system the filter plate is manufactured completely in non-magnetic materials.	
	Filter Cabinet	Cabinet for filter plate inside of the RF cabin for MAGNETOM Symphony/ C!/ Avanto/ Espree/ ESSENZA/ Aera/ Trio/ Verio/ Skyra/ Spectra/ Prisma as well as Biograph mMR. Wall-mounted cabinet custom-made. Outside dimensions	Wall mounted: 

		approx. 1,300 x 420 x 520 mm (L x W x H)	
80.	Emergency hatch	With decorative plate on both sides, white, to be opened from the outside, approx. 600 x 600 mm.	
81.	Observation/ Day-light window	Observation / daylight window that can be opened, installation sideways. Dimensions: 1,200 x 800 mm.	

82.	Threshold-free door	Soundproof threshold-free door for the RF-Cabins of MAGNETOM 1.5T/3T that can be opened and closed without effort by pneumatic mechanism. Clear opening: W x H 1,200 x 2,050 mm	
83.	IV Stand, non-magnetic	Aluminum and plastic, non-magnetic up to 3 tesla. Single-hand height adjustment: 1350 – 2150 mm. 4 plastic bottle hooks: max. 2 kg per hook. Antistatic roller diameter 50 mm, 5 of which lockable	
84.	Instrument table, non-magnetic	Mobile 4-leg frame construction made of 18/10 CrNi steel. Optimal - referring to the handle bar - working height of 800 mm. 63.5 cm x 43.5 cm tubular frame forms a safety rail; no sharp edges. Two shelves: surface ground to reduce glare, resistant to disinfectants, removable. Standard sizes, dimensioned to conform with the DIN standard for sterilized goods. 4 plastic swivel castors, two of them with brakes, wear resistant, non-discoloring when in contact with other material, electrically	

		conductive, Ø 75 mm, 4 wall spacer rings.	
85.	Service Syphon Kit	Service Syphon Kit	 The diagram shows the components of the Service Syphon Kit. It includes a main syphon assembly with a red handle (7) and a black cap (8). There are also several spacer rings (1, 2, 3, 4, 5, 6) and a small assembly (4, 4, 4) shown in a red box. The components are labeled with numbers 1 through 8.
86.	HE-SYPHON DEWAR END	HE-SYPHON DEWAR END	DEWAR SYPHON  The diagram shows the components of the HE-SYPHON DEWAR END. It includes a long, thin syphon tube (1) with a curved end (6). There are also several spacer rings (2, 3, 4, 5) and a small assembly (3) shown in a red box. The components are labeled with numbers 1 through 6.

87.	SYPHON EXT.LEG	SYPHON EXT.LEG	
88.	Quench tube	Quench tube	

89.	Floor plate for table	Floor plate for table	
90.	Adhes. Plastic-Mastic complete	Adhes. Plastic-Mastic complete	
91.	MR-Injector Coupling	MR-Injector Coupling The MR-Injector coupling simplifies synchronized contrast injection and scanning by establishing an active coupling between compatible MR injectors and the MR scanner. The user is guided through an intuitive workflow, with clicks and interactions reduced to the essentials. A direct connection between the MR scanner and the injector control system allows the technologist at the MR operator console to release injections remotely and synchronized with sequences.	 

5.4 Technical Requirements

Technical Requirements based on the key parameters and values as presented in the chapter below.

5.5 Key parameters & values provided by manufacturer

Magnet system

Technical data of the Magnet system:

- Short-bore, patient – friendly design, high homogeneity
- Easy siting due to AS (Active Shielding) and E.I.S. (External Interference Shielding) magnet technology.
- Operating mode: First Level controlled operating mode according to IEC 60601-2-33.

Magnet parameters of MAGNETOM Prisma	
Operating field strength	3 T
Magnet type	Superconductor
Field stability over time	< 0.1 ppm/h
Weight (with cryogenics)	9750 kg
Magnet length	1980 mm ± 2 mm
Open Bore design ¹⁾	600 mm ± 1%
Type of installation	Fixed
Decay characteristics from full field to 20 mT	Approx. 20 s

1) incl. shim coils, gradient coil, RF body coil

Homogeneity (based on highly accurate 24 plane plot)		
	Guaranteed	Typical
10 cm DSV	0.001 ppm	0.0008 ppm
20 cm DSV	0.02 ppm	0.008 ppm
30 cm DSV	0.1 ppm	0.045 ppm
40 cm DSV	0.2 ppm	0.1 ppm
45 cm DSV	0.5 ppm	0.325 ppm
50 cm DSV	1.5 ppm	1.1 ppm

In compliance with the German "Qualifikationsvereinbarung".

Standard deviation V_{rms} (Volume root-mean square) measured with highly accurate 24 plane plot method (20 points per plane).

Standard active shim with 3 linear channels (1st order) and 5 non-linear channels (2nd order).

DSV = Diameter spherical volume (x, y and z direction).

Shimming	
Both: passive and active shimming. Passive shimming during installation.	
Standard active shim with 3 linear channels (1 st order) and 5 non-linear channels (2 nd order).	
3D Shim	Patient-specific automated shim
	Time to shim = approx. 20 s

Shielding		
Active Shielding (AS)	5th generation active shielding (AS) technology with counter coils	
Fringe field (axial×radial)	0.5 mT ¹	6.0 m×3.5 m
	0.1 mT	9.2 m×5.9 m
External Interference Shield (E.I.S.)	Patented shielding system integrated into the magnet	
	Continuous compensation and automatic suppression of external magnetic field interferences during measurement (caused by moving ferromagnetic objects or nearby power lines)	

¹ pacemaker safety limit

Magnet cooling system of MAGNETOM Prisma	
Refill interval (typical) ²	Not applicable
Boil-off rate (typical) ²	0.0 l/year
Max. helium capacity	approx. 933 liters
Minimum helium level	310 l
Cryostat	Stainless steel

² For typical clinical use, depending on sequences and operating time with running helium compressor. The system needs to be serviced at regular interval. Undisturbed magnet cooling for 24 hours and 7 days a week.

Static magnetic field	
Maximum static magnet field	3.1 T
Maximum spatial gradient of static magnet field	7 T/m
Maximum product of the static magnetic field and spatial gradient of static magnet field	17 T ² /m

Patient handling

Fixed table

Max. patient weight for vertical and horizontal table movement	250 kg	
Max. scan range	2050 mm $\pm$ 5 mm	
Vertical table movement	Range	640 mm –1020 mm $\pm$ 13 mm ³
Horizontal table movement	Max. range	2886 mm $\pm$ 10 mm
	Max. speed	200 mm/s $\pm$ 2 mm/s
	Horizontal accuracy for repositioning from one direction	$\pm$ 0.5 mm
Continuous table movement during scan capable		

3 Depending on the floor conditions

Dockable patient table		
Max. patient weight for vertical and horizontal table movement	250 kg	
Max. scan range	2050 mm $\pm$ 5 mm	
Vertical table movement	Range	680 mm –1060 mm $\pm$ 13 mm ³
Horizontal table movement	Max. range	2886 mm $\pm$ 10 mm
	Max. speed	200 mm/s $\pm$ 2 mm/s
	Horizontal accuracy for repositioning from one direction	$\pm$ 0.5 mm
Continuous table movement during scan capable		

3 Depending on the floor conditions

Physiological Measurement Unit (PMU)

Frequency	0.5 Hz - 50 Hz = 30 to 3,000 heartbeats per minute
Modulation parameters	GSFK
Power Options	DC
IP rating	IPX0
Battery characteristics:	
- voltage	3.7 V
- capacity	1320 mAh

Maximim battery life	- Min. 500 charging cycles - Storage: 12 months at room temperature and 50 % SoC lead to minimum recoverable capacity of 85 %
Effective radiated power	- PERU: -3.6 dBm EIRP - PPU: -4.7 dBm EIRP
Frequency/bandwidth of transmission	2402 - 2480 MHz (4 channels, 1 MHz bandwidth, Tx only)
Mode of operation	Continuous operation in accordance with p. 6.6 IEC 60601-1
MR compatibility	MR safe
Type of working part	BF
Dimensions	ECG and respiration measurement module: Base module (electronics): 138 x 57 x 19 mm Heart rate measurement unit: Base module (electronics): 138 x 57 x 19 mm
Weight	ECG sensor / respiratory sensor: 0.265 kg. Heart rate sensor: 0,2 kg.

PERU Sensor Kit:	Battery type	LiPo
	Working time of full-charged battery	12 h – 24 h (depends on the time in bore)
	Time of battery charging	3 h
	Input voltage from charger	5 V (USB)
	Power and current from Charger	500 mA/2.5 W
	Weight	0,610 kg
	Colors of the leads	Grey
	Length of the Leads	Attached cable(s)/hose for connecting ECG electrodes: 500 mm
	MR compatibility	Safe
	Frequency	15 – 190 bpm
	Modulation Type	GFSK
	Parameters of the external trigger signal	Max $\pm$ 5V Min. 5 mA
ECG electrodes:	MR compatibility	MR safe
	Shelf life	24 months from date of manufacture
	Dimensions	Basepad size 5.1 cm square with rounded corners
	Electrode Potential	0.1 mV
	Noise voltage, μ V	4 μ V
	Adhesive parameters	Medical grade acrylic adhesive
	MR compatibility	MR safe
	Special features	Radiotranslucent
	AC Impedance	63 Ω
	Defib Recovery	8.3 mV
	Recovery Rate	0.2 mV/sec
Peripheral Pulse Unit:	Bias	5.2 mV
	Battery type	LiPo
	Working time of full-charged battery	12 h – 24 h (depends on the time in bore)
	Time of battery charging	3 h
	Input voltage from charger	5 V (USB)
	Power and current from Charger	500 mA/2.5 W

	Weight	0,462 kg
	Dimensions	Base unit (electronics): 138 x 57 x 19 mm
	Colors of the leads	Black and grey
	Length of the Leads	Attached flexible cable/tube: length 475 mm
	MR compatibility	safe
	Frequency	15 – 190 bpm
	Modulation Type	GFSK
	Parameters of the external trigger signal	Max $\pm$ 5V Min. 5 mA
PMU charger: (charging station + power supply)	Dimensions	138 x 123 x 35 mm (excl. USB power cable) 55 x 41 x 31,5 mm (power supply)
	Weight	0.604 kg (charging station) 0.390 kg (power supply)
	Power parameters	Max. 5 W
	Input current	Max. 1,000 mA
	Input voltage	100-240V
	Type of voltage	AC
	Frequency	50-60 Hz
	Output voltage	5 V (USB)
	Output current	1400 mA

Hearing protection data

The A-evaluated, effective sound pressure level was measured according to NEMA MS 4-2010(National Electrical Manufacturers Association) using the maximum gradient acoustic noise (MGAN) method.

Patient noise	
XR gradients	Patients require hearing protection with an SNR ⁴ = 14 dB or more

4 SNR: Single number rating

Noise for personnel in the examination room	
XR gradients	Noise measured: 84.4 dB(A)
	Hearing protection: SNR = 2 dB or more

Gradient system

General features	
Actively shielded (AS) whole-body gradient coil system	
Extremely low eddy currents	
Water-cooled coil and amplifier for maximum performance	
All axes force compensated	
Gradient amplifier⁵	
- Water-cooled, highly compact, modular design - Ultra-fast solid-state technology with very low switching losses	
Max. output voltage	2250 V
Max. output current	900 A

Max. power	2.025 MW
------------	----------

5 Values for each of the 3 gradient axes

Gradient performance for each axis	
Max. amplitude	80 mT/m
Min. rise time	400 μ s
Max. slew rate	200 T/m/s
Vector gradient performance (vector addition of all 3 gradient axes)	
Max. eff. amplitude	139 mT/m
Max. eff. slew rate	346 T/m/s
Gradient duty cycle	100%

RF technology

Transmit technology		
Frequency stability (5 min)	$\pm 2 \times 10^{-10}$	
Frequency control	32 bits (0.03 Hz)	
Phase control	16 bits (0.006 degrees)	
Body coil	Integrated whole body no tune transmit/receive coil with 32 rungs	
	Optimized RF efficiency and signal-to-noise ratio (SNR)	
	Real time feedback loop for unmatched RF stability	
Transmitter path	Feedback loop for excellent RF stabilization	
	Transmit amplitude	16 bit control 25 ns resolution
	Gain stability (after first minute)	<0.05 dB (1 s) <0.2 dB (5 min)
Transmit amplifier	Extremely compact, water-cooled solid state amplifier, fully integrated at the magnet as part of DirectRF™ technology	
	Transmit amplifier bandwidth	500 kHz
Peak power	(TimTX TrueShape 2ch transmit)	
Power capability of the amplifier	43.2 kW	
Channel 0	17.5 KW	
Channel 1	25.7 KW	
TimTX TrueForm		
• Innovative techniques in the RF excitation hardware • New application and processing features enabling uniform RF distribution in all body regions		
TrueForm excitation uses amplitude and phase transmission settings optimized for dedicated body regions. Feeding two independent ports of the integrated bodycoil with an optimized weighting yields a homogeneous B1 distribution		
B1 Filter	An adaptive inline image filter that reduces any remnant B1 effects without affecting image contrast	

RF receiver technology	
Maximum number of channels ⁶	204
Number of independent receiver channels	48 ⁷ , 64, 128 ⁸
Quadrature demodulation and filtering	Digital
Receiver bandwidth	500 Hz–1 MHz (for each channel)
Receiver signal resolution	32 bit
ADC sampling rate	80 MHz
Preamplifier noise figure	< 0.5 dB
Dynamic range at coil connector (referred to 1 Hz resolution bandwidth)	164 dB instantaneous at receiver 169 dB with automatic gain control at local coil connector

6 Channels (coil elements) that can be connected simultaneously

7 Only available for MAGNETOM Prisma^{fit}

8 Optional

Multi-Nuclear Option ⁹				
- Multi-Nuclear Imaging and Spectroscopy for the following nuclei: ³He, ⁷Li, ¹³C, ¹⁷O, ¹⁹F, ²³Na, ³¹P, ¹²⁹Xe - Decoupling nuclei: ¹³C, ³¹P - Hyperpolarization nuclei: ³He, ¹³C, ¹²⁹Xe				
Frequency band	³ He	⁷ Li	¹³ C	¹⁷ O
Center frequency [MHz]	93.8524	47.8799	30.9806	16.7012
Usable transmitter band width [kHz]	±50	±50	±50	±50
Usable receiver band width [kHz]	±190	±100	±65	±50
	¹⁹ F	²³ Na	³¹ P	¹²⁹ Xe
Center frequency [MHz]	115.9040	32.5885	49.8719	34.0816
Usable transmitter band width [kHz]	±50	±50	±50	±50
Usable receiver band width [kHz]	±235	±65	±100	±70

9 Optional

RF coils

Coil name	Coil type	Coil design	Intended use		General data	
			Resonant nuclei	Applications	Dimensions	Weight
Body 18 Body 18 long	Receive	18 channel phased array	¹ H	- Head - Neck - Thorax - Heart - Abdomen - Pelvis - Hip - Vascular	385 x 590 x 65 mm	1,6 kg (in configuration with a long cable 2 kg)

Coil name	Coil type	Coil design	Intended use		General data	
			Resonant nuclei	Applications	Dimensions	Weight
Body 30	Receive	30 channel phased array	^1H	• Thorax • Heart • Abdomen • Pelvis • Hip • Vascular	300 x 390 x 315 mm	3 kg
Head 32	Receive	32 channel phased array	^1H	• Head examination	300 x 390 x 315 mm	5.6 kg
Head/Neck 20	Receive	20 channel phased array	^1H	• Headexamination • Neckexamination • MRHead/Neck Angiography • Combined head/neck examination • TMJ (temporo-mandibular joints) • Proton spectroscopy	440 x 330 x 370 mm	4.7 kg
Head/Neck 64	Receive	64 channel phased array	^1H	• Headexamination • Neckexamination • MRHead/Neck Angiography • Combined head/neck examination • TMJ (temporo-mandibularjoints) • Proton spectroscopy	435 x 395 x 350 mm	5.2 kg
Spine 32	Receive	32 channel phased array	^1H	• Spine imaging • Imaging upper body, abdomen, pelvis, liver, and heart (when combined with an anterior coil)	1200 x 489 x 63 mm	11 kg
Flex Large 4	Receive	4 channel phased array	^1H	• Imaging of large regions such as medium to large shoulder, hip, and knee	516 x 224 mm	0.55 kg
Flex Small 4	Receive	4 channel phased array	^1H	• Imaging of small regions such as small to medium shoulder, wrist, elbow, and ankle	366 x 174 mm	0.45 kg
Peripheral Angio 36	Receive	36 channel phased array	^1H	• High resolution angiography of both legs with high signal-to-noise ratio • Bilateral examinations of long bones of the legs	860 x (300 up to 640) x 280 mm	8 kg
Loop 11	Receive	Linear	^1H	• Examination of upper or lower extremities (e.g. shoulder, axilla)	Diameter: 110 mm	0.225 kg
Loop 7	Receive	Linear	^1H	• Examination of inner ear, structure of wrist and fingers, pediatrics examinations¹⁰	Diameter: 70 mm	0.2 kg
Loop 4	Receive	Linear	^1H	• Examination of small structures near the	Diameter: 40 mm	0.175 kg

Coil name	Coil type	Coil design	Intended use		General data	
			Resonant nuclei	Applications	Dimensions	Weight
				surface, e.g. joints of fingers and toes, wrist, skin, temporomandibular joints (TMJ)		
Hand/Wrist 16	Receive	16 channel phased array	¹ H	High resolution hand and wrist imaging	Coil: 332 × 215 × 115 mm Base plate: 527 × 470 × 55 mm	Coil: 2.8 kg Base plate: 1.6 kg
Foot/Ankle 16	Receive	16 channel phased array	¹ H	High resolution foot and ankle imaging	Coil: 410 × 330 × 390 mm Base plate: 427 × 333 × 383 mm	Coil 3.2 kg Base plate: 7.1 kg
Shoulder Large 16 Shoulder Small 16	Receive	16 channel phased array	¹ H	High resolution shoulder imaging	Hole Large 200 mm Hole Small 165 mm Base plate: 445 × 490 mm	Large: 2.3 kg Small: 2.2 kg Base plate: 1.7 kg
UltraFlex Large 18	Receive	18 channel phased array	¹ H	- Hand and wrist - Elbow and shoulder - Hip and knee - Foot and ankle - Head and neck	590 × 290 × 14 mm	1.8 kg
UltraFlex Small 18	Receive	18 channel phased array	¹ H	- Hand and wrist - Elbow and shoulder - Hip and knee - Foot and ankle - Head and neck	410 × 190 × 14 mm	1.4 kg
TxRx Knee 15	Transmit/Receive	Transmit: quad-rature Receive: 15 Channel phased array	¹ H	- Examinations of joints in the area of the lower extremities - High resolution knee imaging	500 × 537 × 305 mm Minimum inside diameter: 154 mm	6.2 kg
TxRx Knee 15 Flare	Transmit/Receive	Transmit: quad-rature Receive: 15 Channel phased array	¹ H	- Examinations of joints in the area of the lower extremities - High resolution knee imaging	Height (coil + base plate): 305 mm Coil: 355 × 270 mm Base plate: 500 × 537 mm	6.2 kg (including base plate)
TxRx 15 Ch Knee	Transmit/Receive	Transmit: quad-rature Receive: 15 Channel phased array	¹ H	- Examinations of joints in the area of the lower extremities - High resolution knee imaging	500 × 537 × 305 mm Minimum inside diameter: 154 mm	6.2 kg
4 Ch BI Breast	Receive	4 Channel phased array	¹ H	Simultaneous basic imaging of both breasts in all directions	880 × 470 × 210 mm	10 kg

Coil name	Coil type	Coil design	Intended use		General data	
			Resonant nuclei	Applications	Dimensions	Weight
				• Uni-lateral biopsy imaging for lateral and medial access • High resolution 2D and 3D imaging • For quantitative spectroscopy (syngo GRACE) a reference bottle can be inserted		
16 Ch BI Breast	Receive	16 Channel phased array	¹ H	• Simultaneous basic imaging of both breasts in all directions • Uni-lateral biopsy imaging for lateral and medial access • High resolution 2D and 3D imaging • For quantitative spectroscopy (syngo GRACE) a reference bottle can be inserted	Frame: 880 x 470 x 210 mm Coil: 550 x 316 x 166 mm	Coil: 5.8 kg Including frame: 14 kg
Breast 18	Receive	18 Channel phased array	¹ H	• Simultaneous imaging of both breasts in all directions • Axillary imaging elements • High-resolution 2D and 3D imaging • For quantitative spectroscopy (syngo GRACE) a reference bottle can be inserted	575 x 410 x 205 mm	5.5 kg
2-/4-/8-Channel Sentinel Breast Coil	Receive	Optional 2, 4, or 8 Channel phased array	¹ H	• Simultaneous imaging of both breasts in all directions • Uni- or bi-lateral imaging of the breasts in sagittal direction • Uni-lateral biopsy imaging for lateral and medial access • High-resolution 2D and 3D imaging • For quantitative spectroscopy (syngo GRACE) a reference bottle can be inserted	1097 x 582 x 279 mm	16 kg With plate: 22 kg

Coil name	Coil type	Coil design	Intended use		General data	
			Resonant nuclei	Applications	Dimensions	Weight
2-/10-/16-Channel Sentinel Breast Coil	Receive	Optional 2, 10, or 16 Channel phased array	¹ H	• Simultaneous imaging of both breasts in all directions • Uni- or bi-lateral imaging of the breasts in sagittal direction • Uni-lateral biopsy imaging for lateral and medial access • High-resolution 2D and 3D imaging • For quantitative spectroscopy (syngo GRACE) a reference bottle can be inserted	1097 x 582 x 279 mm	16 kg With plate: 22 kg
Endorectal Coil Kit (for connecting the disposable endorectal coil; to be ordered separately)	Receive	Linear	¹ H	• Prostate • Colon • Cervix Uteri	Adapter: 690 mm x 45 mm x 35 mm Cushion: 370 mm x 301 mm x 175 mm	0.5 Kg 1.9Kg
TxRx CP Head	Transmit/Receive	Quadrature	¹ H	• Head examinations	315 x 475 x 360 mm	5.2 kg
TxRx CP Extremity	Transmit/Receive	Linear	¹ H	• Knee • Ankle	256 x 360 x 310 mm	6.6 kg
Special Purpose 4	Receive	4 Channel phased array	¹ H	• Carotids • Examinations with small Field-of-View • Small structures near the surface	143 x 90 x 20 mm	0.2 kg
Pediatric 16	Receive	16 Channel phased array	¹ H	• Imaging of neonates and infants ¹⁰⁾	Coil: 371 x 318 x 269 mm Cradle: 685 x 319 x 123 mm	Coil: 4.4 kg Front coil: 1.4 kg Cradle: 2.3 kg

¹⁰ MR scanning has not been established as safe for imaging fetuses and infants under two years of age. The responsible physician must evaluate the benefit of the MRI examination in comparison to other imaging procedures.

Sequences

Spin echo family of sequences

- SE (Spin Echo): single and multi echo
- TSE (Turbo Spin Echo)
- HASTE (Half-Fourier Acquisition with Single Shot Turbo Spin Echo)
- SPACE (Sampling Perfection with Application-optimized Contrast using different flip angle Evolutions): for fast 3D spin echo imaging

Gradient echo family of sequences

- FLASH (Fast Low Angle SHot): spoiled GRE (gradient echo)
- VIBE (Volume Interpolated Breathhold Examination): for fast 3D gradient echo imaging
- MEDIC (Multi Echo Data Image Combination)
- DESS (Dual Echo Steady State)
- TurboFLASH
- FISP (Fast Imaging with Steady State Precession)
- EPI (Echo Planar Imaging): with signal type SE or FID (Free Induction Decay)
- Segmented EPI: segmented in phase or readout direction
- ToF (Time-of-Flight) Angiography
- PC (Phase Contrast) Angiography
- CV/BEAT: comprehensive sequence for cardiac imaging

Spectroscopy sequences

- SVS (Single Voxel Spectroscopy)
- FID (Free Induction Decay)
- CSI (Chemical Shift imaging)

Computer system

The specifications contained in this section apply to new system deliveries with XA30A and depending on appropriate hardware preconditions also to systems which are upgraded to XA30A.

Host computer	Processor	Intel Xeon $\geq$ W-2133 (6 Core)
	Clock rate	3.6 GHz or comparable
	Main memory (RAM)	64 GB
	Hard disk (DICOM Standard, ISO 9660)	SSD: 480 GB Approx. Images: 256 ¹¹⁾
	CD/DVD drive	Not built in, but optionally connectable by USB
	Media drive	SDHC card reader
Advanced host computer	Processor	Intel Xeon $\geq$ W-2145 (8 Core)
	Clock rate	3.7 GHz or comparable
	Main memory (RAM)	96 GB
	Hard disk (DICOM Standard, ISO 9660)	SSD: 480 GB Approx. Images: 256 ¹¹⁾
	CD/DVD drive	Not built in, but optionally connectable by USB
	Media drive	SDHC card reader
Color LCD monitor ¹²⁾ High resolution widescreen monitor:	Screen size (diagonal)	24"
	Horizontal scanning frequency	31–76 kHz
	Vertical scanning frequency	59–61 Hz
	Screen matrix	1920×1200 pixels

syngo MR Workplace	• Horizontally tiltable, forward and backward • Automatic backlight control for long-term brightness stability	
	Processor	Intel Pentium J5005 (4 Core)
	Clock rate	1.5 GHz
	Main memory (RAM)	8 GB
	Hard disk	SSD: 128 GB
	CD/DVD drive	Not built in, but optionally connectable by USB
	Audio	On board
	Ethernet	On board, 1 Gbit

- 11 Using Enhanced DICOM > 2 200 000 image stacks with 25 slices and a matrix size of 256×256 can be stored
- 12 A standard monitor without calibration is not suitable for diagnostic purposes. Please consider the initial acceptance testing for image display devices and the follow-up service for constancy testing on a regular base, as offered by Siemens service.

Measurement and reconstruction system for Tim [204 × 64]	Processor	2 x Intel Xeon E5-2620v4 (8 Core)
	Clock rate	2 x 2.1 GHz or comparable
	Main memory (RAM)	64 GB
	SSD for raw data	≥ 480 GB
	SSD for system software	≥ 240 GB
	Reconstruction speed	40 404 recons per second (256 ² FFT, full FoV) 149 532 recons per second (256 ² FFT, 25% recFoV)
	Parallel scan and recon	Simultaneous scan and reconstruction of up to 12 data sets
	GPGPU	1 x Nvidia Quadro P4000
Advanced Measurement and reconstruction system for Tim [204 × 64] Option 1	Processor	2 x Intel Xeon E5-2640v4 (10 Core)
	Clock rate	2 x 2.4 GHz, or comparable
	Main memory (RAM)	128 GB
	SSD for raw data	≥ 2 x ≥ 480 GB
	SSD for system software	≥ 240 GB
	Reconstruction speed	52 493 recons per second (256 ² FFT, full FoV) 224 719 recons per second (256 ² FFT, 25% recFoV)
	Parallel scan and recon	Simultaneous scan and reconstruction of up to 12 data sets
	GPGPU	2 x Nvidia Quadro P4000
	Processor	2 x Intel Xeon E5-2680v4 (14 Core)
	Clock rate	2 x 2.4 GHz, or comparable
	Main memory (RAM)	256 GB

Advanced Measurement and reconstruction system for Tim [204 × 64] Option 2	SSD for raw data	2 × ≥ 480 GB
	SSD for system software	≥ 240 GB
	Reconstruction speed	77 994 recons per second (256 ² FFT, full FoV) 187 919 recons per second (256 ² FFT, 25% recFoV)
	Parallel scan and recon	Simultaneous scan and reconstruction of up to 12 data sets
	GPGPU	2 × Nvidia Quadro P4000
Measurement and reconstruction system for Tim [204 × 128]	Processor	2 × Intel Xeon E5-2640v4 (10 Core)
	Clock rate	2 × 2.4 GHz, or comparable
	Main memory (RAM)	128 GB
	SSD for raw data	≥ 2 × ≥ 480 GB
	SSD for system software	≥ 240 GB
	Reconstruction speed	52 493 recons per second (256 ² FFT, full FoV) 224 719 recons per second (256 ² FFT, 25% recFoV)
	Parallel scan and recon	Simultaneous scan and reconstruction of up to 12 data sets
Advanced measurement and reconstruction system for Tim [204 × 128]	GPGPU	2 × Nvidia Quadro P4000
	Processor	2 × Intel Xeon E5-2680v4 (14 Core)
	Clock rate	2 × 2.4 GHz, or comparable
	Main memory (RAM)	256 GB
	SSD for raw data	2 × ≥ 480 GB
	SSD for system software	≥ 240 GB
	Reconstruction speed	77 994 recons per second (256 ² FFT, full FoV) 187 919 recons per second (256 ² FFT, 25% recFoV)
	Parallel scan and recon	Simultaneous scan and reconstruction of up to 12 data sets
	GPGPU	2 × Nvidia Quadro P4000

Euro Connector (System IEC 320)	
Load rating	100–240 Vac
	50–60 Hz
	max. 10 A

Network/data coupling

Network/data coupling

The MR system provides the level of safety according IEC 60950-1 outside patient environment.
All equipment connected to the system's network/data couplings must also provide minimum level of safety according IEC 60950-1.

Basic syngo MR viewing, post-processing, and filming functionalities

Viewing

- Image display
 - Various display layouts selectable
 - Up to 4 patients can be simultaneously active
 - Image annotation and labeling
 - Configurable image text
 - Interactive movie
- Windowing
 - Freely selectable window width and center
 - Windowing (on all, succeeding, or single images)
 - Auto-windowing for optimized contrast ROI/VOI evaluation
- ROI/VOI evaluation

Parallel evaluation of multiple regions/volumes of interest

 - Circle
 - Freehand
 - Statistical evaluation:
 - Area/Volume
 - Standard deviation
 - Mean value
 - Min/max values
- Range creation
 - New reconstructions based on acquired images
 - Parallel, radial, radial sliced, and curved ranges
- Position display
 - Displays measured slice positions on localizer image and selected series
- Image manipulations
 - Reversal of gray-scale values
 - Image rotation by user-defined angle
 - Image zoom and pan

Mean Curve

- Time-intensity analysis

Calculations

Arithmetic operations on images and series:

- Addition, subtraction, multiplication, division within and across series
- Calculated b-value and ADC maps
- Image filter for smoothing and sharpening images
- Motion correction (elastic) between volumes or phases
- 2D and 3D distortion correction

Composing/Combine

- Composing of sagittal or coronal images (for example, whole-spine exams or run-offs)
- Combine different image stacks into a single one

3D

- MPR (Multi-Planar Reconstruction): Real-time multi-planar reformatting of secondary views
- MIP (Maximum Intensity Projection): 3D reconstructions of vessels from a 3D data set, or a 2D sequential slice data set (acquired with dedicated MR Angiography sequences)
- MinIP (Minimum Intensity Projection): Similar to MIP but reconstructs the minimum intensity (e.g. for Dark Blood techniques)
- VRT (Volume Rendering Technique): 3D rendering with specific adaptable opacity and shading

Filming

- Connection via DICOM Basic Print
- Selectable positioning of images onto virtual film sheet
- Selectable various film layouts
- Windowing, image zoom and pan on film sheet

Specific absorption rates

The specific absorption rate for the RF power and the rate of change of the gradient fields are checked according to the requirements of norm IEC 60601-2-33.

Environmental conditions

Control room	
Temperature	15–30 °C
Relative air humidity	40 to 60%
Absolute air humidity	<11.0 g/kg
Air pressure	1060 hPa–700 hPa

Electronics room	
Temperature	15–30 °C
Relative air humidity	40 to 80%
Absolute air humidity	<11.0 g/kg
Air pressure	1060 hPa–700 hPa

Examination room	
Temperature	18–22 °C
Relative air humidity	40 to 60%
Absolute air humidity	<11.0 g/kg
Air pressure	1060 hPa–700 hPa

Specific absorption rates

The specific absorption rate for the RF power and the rate of change of the gradient fields are checked according to the requirements of norm IEC 60601-2-33.

Cooling system

Two different customer specific cooling alternatives (Separator or Eco Chiller) are available.

Separator option (for connection to available cooling system)	Water consumption	100 l/min ± 10 l/min ¹¹
	Heat dissipation to water	60 kW
Eco Chiller option with automatic adaptation to the required cooling demands (e.g. different night/day mode) to decrease energy cost ¹⁴	GREEN Cooling Package¹²: - Automatic start if the surrounding temperature is 18 °C (64 °F) or less - If the temperature is less than –10 °C (14 °F) the chiller is switched off¹³	

11 Water temperature: 6°C–14°C (43°F–57°F); allowed delta T: ± 2 °C with max. 1 °C / 30 s

12 Free Cooling Unit, optional

13 In case of clinical routine measurement conditions

14 only available for XR gradients

Line power supply

XR gradients	Values	Tolerance
Voltage	380 V, 400 V, 420 V, 440 V, 460 V, 480 V	±10%
Frequency	50/60 Hz	±1 Hz
Connection value	93 kVA	

For all products, line power has to be supplied via an on-site system contact or via another multipole shut-down mechanism. Room installation has to be in compliance with VDE 0100-710.

In all countries, compliance with local and national legal regulations is required. However, we strongly recommend compliance with the regulations described herein – to the extent permitted by relevant local and national laws – in order to ensure the safety of operating personnel, patients, and third parties.

Power rating plates

XR gradients

380/400/420/440/460/480	V ~
3PHASE	50/60 Hz
SHORT – TIME	kVA: 87
LONG – TIME	kVA: 29

Power consumption

XR gradients

System off ¹⁴	6.1 kW
System ready to measure ¹⁴	8.7 kW
Scan ¹⁴	18.9 kW
Remarks:	
• All values are typical values, applicable for 400 V/50 Hz. • Consumption for optional separator pump not included.	

- 14 The power consumption described herein is based on results that were achieved in a setting according to the COCIR methodology MRI Measurement of the energy consumption (<http://www.cocir.org/index.php?id=46>). Since many variables impact power consumption (e.g. sequences used for scanning and sequence parameters, scan time), there can be no guarantee that each customer will achieve the same values

Dimensions

Component	Width [m]	Depth [m]	Height [m]	Weight [kg]
Examination				
Magnet 3 Tesla AS (incl. Helium)	2.16	1.98	2.54	9750
Magnet in operation, incl. gradient coil, body coil, Tim table, and covers	2.55	4.72	2.54	12250
Tim table	0.76	2.49	0.52 –1.02; +0.013 ¹⁵	
Required min. room height clearance			2.4 ¹⁶	
Min. transport dimensions	2.40	3.58	2.42	
Control Room				
syngo Acquisition Workplace (table+monitor)	1.2	0.8	1.17 (0.72+0.45)	
Host computer	0.22	0.46	0.47	
syngo MR Workplace (optional) (table+monitor)	1.2	0.8	1.17 (0.72+0.45)	
Equipment Room				
Electronics cabinet, incl. system control, RF system, gradient power system, image processor	1.6	0.65	1.98 ¹⁷	1500
Cooling system	0.65	0.65	1.89	500
Heat dissipation	≤5 kW, only ventilation might be required			

15 Depending on the floor conditions

16 Finished floor to finished ceiling

17 Without attachments

Image parameters

Image Parameters	Characteristics
Signal to noise ratio (SNR)	Not worse than 175
High contrast resolution (2D)	Not more than 0.8 mm
High contrast resolution (3D)	Not more than 0.8 mm
Slice thickness	5 ± 1 mm
Low contrast resolution	Not worse than 3 mm

Software

Name and Version	syngo MR XA30
Release date	07.04.2021


5.6 Labeling information

At the electronics cabinet

Product identification label of MAGNETOM Prisma:
MDD (Medical Device Directive):

MDR (Medical Device Regulation):

(1) Configurable device

A configurable device is a device that consists of several components which can be assembled by the manufacturer in multiple configurations.


Those individual components may be devices in themselves.

(2) Identification as medical device

(3) Serial number of the system

(4) Unique device identification

Combination label (MDR)

This label indicates that there are approved combinations of the MR system described in this manual with other products. Please contact Siemens Healthcare GmbH if you have any questions regarding permitted combinations.

Labels according to EN 50419:2006

Within the EU, products identified with this symbol are subject to guidelines 2002/96/EC for old electrical or electronics system, modified by guidelines 2012/19/EU respectively. Please contact Siemens Service in case of questions about returning and disposing the MR system and/or its components and accessories.

Revision label

Revision:														Model - No.:																		
X	X	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65
66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98


The current revision is marked with an X in the respective field. The example above shows revision level 01 of the respective component.

On the magnet cover

Product identification label of MAGNETOM Prisma:
MDD (Medical Device Directive):

MDR (Medical Device Regulation):

(1) Configurable device

A configurable device is a device that consists of several components which can be assembled by the manufacturer in multiple configurations.


Those individual components may be devices in themselves.

(2) Identification as medical device

(3) Serial number of the system

(4) Unique device identification

Combination label (MDR)

This label indicates that there are approved combinations of the MR system described in this manual with other products. Please contact Siemens Healthcare GmbH if you have any questions regarding permitted combinations.

At the magnet bore (front)

Risk of injury

Observe operator manual

6.10.1 For laser light localizer

Safety certificate label (International)

EAC Label*

Each unit of the MAGNETOM Prisma delivered to Russia, Belarus and Kazakhstan will be endorsed with the EAC label right beside the system labels.

**This symbol is applicable to the system until 12/11/2022, after this time the symbol is not applied.*

Label in Russian language for MAGNETOM Prisma:

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
On the coils and interfaces

Name of component	Labels
Head/Neck 20 coil	 The label for the Head/Neck 20 MR Coil 3T includes the following information: SIEMENS Head/Neck 20 MR Coil 3T - GTIN: (01) 04056869014869 - MODEL: (240) 10496500 - ORIGIN DE: (422) 276 - SERIAL: (21) 30744 - IVK - Siemens Healthcare GmbH, Henkestr. 127, 91052 Erlangen, Germany - Made in Germany - ES: 06 01 02 03 04 05, 06 07 08 09 10 11 12, 13 14 15 16 17 18 19 152003 - Safety symbols: Person with lightning bolt, Warning triangle.

Spine 32 coil	 The label for the Spine 32 MR Coil 3T features the Siemens logo at the top. Below it, the product name 'Spine 32 MR Coil 3T' is printed. A QR code is located on the left side. To the right of the QR code, the following information is listed: (01) 04056869014837 GTIN, (240) 10496545 MODEL, (422) 276 ORIGIN DE, and (21) 30688 SERIAL. At the bottom, it states 'Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany' and 'Made in Germany'. A small icon of a person is visible in the bottom left corner of the label area.
Body 18 Body 18 long	Body 18  The label for the Body 18 MR Coil 3T features the Siemens logo at the top. Below it, the product name 'Body 18 MR Coil 3T' is printed. A QR code is located on the left side. To the right of the QR code, the following information is listed: (01) 04056869014876 GTIN, (240) 10496515 MODEL, (422) 276 ORIGIN DE, and (21) 30237 SERIAL. At the bottom, it states 'Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany' and 'Made in Germany'. A small icon of a person is visible in the bottom left corner of the label area. Body 18 long  The label for the Body 18 long MR coil 3T features the Siemens logo at the top. Below it, the product name 'Body 18 long MR coil 3T' is printed. A QR code is located on the left side. To the right of the QR code, the following information is listed: (01) 04056869015392 GTIN, (240) 10835253 MODEL, (422) 276 ORIGIN DE, and (21) 30107 SERIAL. At the bottom, it states 'Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany' and 'Made in Germany'. A small icon of a person is visible in the bottom left corner of the label area.

Flex Small/Large 4	Flex Large 4  SIEMENS Flex Large 4 MR coil 3T (01) 04056869015439 GTIN (240) 10835335 MODEL (422) 276 ORIGIN DE (21) 30196 SERIAL IVK Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany Made in Germany  Flex Small 4  SIEMENS Flex Small 4 MR coil 3T (01) 04056869014395 GTIN (240) 08625779 MODEL (422) 276 ORIGIN DE (21) 30241 SERIAL IVK Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany Made in Germany 
UltraFlex small/large 18	UltraFlex large 18  SIEMENS Flex Large 18 MR Coil 3T (01) 04056869044705 GTIN (240) 10835680 MODEL (422) 276 ORIGIN DE (21) 1006 SERIAL IVK Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany Made in Germany 

	UltraFlex small 18 
Head/Neck 64	
Tx/Rx CP Head Coil	


Body 30/60	 SIEMENS Body 30 MR coil 3T (01) 04056869015460 GTIN (240) 10835446 MODEL (422) 276 ORIGIN DE (21) 30120 SERIAL IVK Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany Made in Germany
Peripheral Angio 36	 SIEMENS PA 36 MR Coil 3T (01) 04056869014913 GTIN (240) 10496543 MODEL (422) 276 ORIGIN DE (21) 30131 SERIAL IVK Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany Made in Germany SIEMENS (1P) Model No. 10496543 (S) Serial No. 1081 IVK Siemens AG, Wittelsbacherplatz 2, DE-80333 Muenchen Made in Germany 612988 613155 ES: 06 07 08 09 10 11 12 13 14 15 16 17 18 19
TxRx 15 Ch Knee TxRx Knee 15 TxRx Knee 15 Flare	TxRx 15 Ch Knee



TxRx Knee 15

TxRx Knee 15 Flare

Label Part No. 6000748

Tx/Rx CP Extremity Coil	
Hand/Wrist 16	 
Foot/Ankle 16	 

Head 32

Pediatric 16

 Quality Electro Dynamics, LLC
6655 Beta Drive, Suite 100, Mayfield Village, OH 44143 USA
Distributed by Siemens Healthcare GmbH

MD

Pediatric 16, A 3T TIM Coil (Receive Only)

IVK

Siemens	REF	(1P) 10835523
QED	REF	Q7000137
SN	(S)	NNNN
Freq. MHz		123.2
		YYYY-MM-DD

ES	00	02
03	04	05
06	07	08
09	10	11

CE 2797

EC REP

EMERGO EUROPE
Pruissengracht 20
2524 AP The Hague
The Netherlands

RX ONLY

(01)00814848020577
(11)YYMMDD
(21)NNNN
(240)10835523
(422)B40

Pot. www.QEDPN.com

Label No. 6000036

Multi-nuclear spectroscopy coil

Head Coil $^{13}\text{C}^1\text{H}$

Head Coil $^{31}\text{P}^1\text{H}$

Head Coil $^{23}\text{Na}^1\text{H}$

Manufactured by
RAPID Biomedical GmbH
Kettelerstr. 3-11, 97222 Rimpf
Germany - www.rapidbiomed.co

$^{31}\text{P} / ^1\text{H}$
Head Coil 3T

Model Type	V-XQ-HQ-030	
REF	01487	
SN	027	
Version	V02	
	2019-11	
Nucleus	^1P	^1H
Frequency	49.9 MHz	123 MHz
Channels	TxRx (quadrature)	TxRx (quadrature)
MR System	Siemens MAGNETOM 3T	

CE 0197

ATTENTION: Do not use if there is damage of enclosure, connections or cord.

Coils for mammary glands with a modular design - no more than 2 pcs. (if necessary), versions:

- 8-ch. Sentinelle Breast Coil;
- and/or
- 16-ch Sentinelle Breast Coil.

16Ch

16 CH BREAST SEN MR COIL KIT 3T

Invivo Corporation
3545 SW 47th Avenue
Gainesville, FL 32605-7891
U.S.A. +1 (352) 336-0910

REF 
(1P) 10 60 66 10

SN 
(S) serial

Rev 
(2P) 03

Distributed by:
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany

EC REP

0123

GBM Authorized
Representative Ltd
The White House
2 Meadow
Godalming GU7 3HN
United Kingdom
+44 (0)7710 039 721
www.mba-gbm.com

8Ch

8CH BREAST SEN MR COIL KIT 3T

Invivo Corporation
3545 SW 47th Avenue
Gainesville, FL 32605-7891
U.S.A. +1 (352) 336-0910

REF 
(1P) 10 60 66 03

SN 
(S) serial

Rev 
(2P) 05

Distributed by:
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany

0123

Coils for mammary glands with a modular design for imaging and/or biopsy - not more than 2 pcs. (if necessary), versions:
- 2-/4-/8-Channel Sentinelle Breast Coil;
and/or
- 2-/10-/16-Channel Sentinelle Breast Coil

2-4 Ch

Sentinelle
by Invivo

2/4CH BREAST SEN MR COIL KIT 3T

REF 
(1P) 10 60 66 04

SN 
(S) serial

Rev 
(2P) 04

Invivo Corporation
3545 SW 47th Avenue
Gainesville, FL 32608-7891
U.S.A. +1 (352) 336-5810

YYMMDD

Distributed by:
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany

CE **EC** **REP**
0123

GBM Authorized
Representative Ltd
The White House
2 Meadow
Godalming GU7 3HN
United Kingdom
+44 (0)7710 039 721
www.mba-gbm.com

16Ch

Sentinelle
by Invivo

16 CH BREAST SEN MR COIL KIT 3T

REF 
(1P) 10 60 66 10

SN 
(S) serial

Rev 
(2P) 03

Invivo Corporation
3545 SW 47th Avenue
Gainesville, FL 32608-7891
U.S.A. +1 (352) 336-5810

YYMMDD

Distributed by:
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany

CE **EC** **REP**
0123

GBM Authorized
Representative Ltd
The White House
2 Meadow
Godalming GU7 3HN
United Kingdom
+44 (0)7710 039 721
www.mba-gbm.com

8Ch

	Invivo Corporation 3545 SW 47th Avenue Gainesville, FL 32608-7691 U.S.A. +1 (352) 336-0910 8CH BREAST SEN MR COIL KIT 3T REF (1P) 10 60 66 03 SN (S) serial Rev (2P) 05 Distributed by: Siemens Healthcare GmbH Henkestr. 127 91052 Erlangen Germany CE 0123
Breast 18	SIEMENS Breast 18 MR coil 3T (01) 04056869015286 GTIN (240) 10606916 MODEL (422) 276 ORIGIN DE (21) 30117 SERIAL IVK Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany Made in Germany
4-channel BI Breast Coil 16Ch AI Breast coil	4 CH BI Breast MR Coil Kit 3T Medical equipment for MRI imaging Manufacturer: KORAS MRI products GmbH Lehrstraße 4 97204 Hirschberg Germany Distributed by: Siemens Healthcare GmbH Koras Ref.: 113113 Rev. 06 (01) 040512800508 (422) 276 (21) 30117 (040) 1015528 GTIN Origin Rev Serial Model Notes (label for repeat order) CE
	16-Channel AI Breast Coil

	16Ch AI Breast Coil 3T RAPID Biomedica Manufactured by: RAPID Biomedica GmbH Kettelerstr. 3-11, 91222 Regensburg, Germany www.rapidbiomed.de Distributed by: Siemens Healthcare GmbH Herkenstraße 127 91052 Erlangen Germany  (01) 04260487680041 (20) 02 (21) 1433 (240) 10185782
FlexCoil Interface 3T	Flex coil interface 3T  (01) 04056869154015 GTIN (240) 10500091 MODEL (422) 156 ORIGIN CN (21) 4123 SERIAL IVK Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center, Gaoxin C. Ave., 2nd, Hi-Tech Industrial Park, 518057 Shenzhen, PEOPLE'S REPUBLIC OF CHINA Made in China 
Tim Coil Interface	TIM Coil Interface 3.0T  (240) 10500092 MODEL (422) 156 ORIGIN CN (21) 4651 SERIAL IVK 
Tim Coil Interface MNO	MNO Coil Interface 3.0T  (240) 10656124 MODEL (422) 156 ORIGIN CN (21) 1392 SERIAL IVK 


Multi-channel Interface #3T	SIEMENS Multi – Channel Interface 3T (01) 04056869119162 GTIN (240) 11134266 MODEL (422) 276 ORIGIN DE (21) 1001 SERIAL IVK Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany Made in Germany 
Loop Coil, large	SIEMENS 3T Loop 11cm MR coil (01) 04056869014357 GTIN (240) 08625662 MODEL (422) 276 ORIGIN DE (21) 30037 SERIAL IVK Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany Made in Germany 
Loop Coil, medium	SIEMENS 3T Loop 7 cm MR coil (01) 04056869014364 GTIN (240) 08625670 MODEL (422) 276 ORIGIN DE (21) 30052 SERIAL IVK Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany Made in Germany 

Loop Coil, small	 SIEMENS 3T Loop 4 cm MR coil (01) 04056869014371 GTIN (240) 08625720 MODEL (422) 276 ORIGIN DE (21) 30055 SERIAL Siemens Healthineers GmbH, Henkestr. 127 91052 Erlangen, Germany Made in Germany IVK 
Special Purpose Coil	 special purpose coil 4 (01) 04056869015330 GTIN (240) 10654891 MODEL (422) 156 ORIGIN CN (21) 1503 SERIAL Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center, Gaoxin C. Ave., 2nd, Hi-Tech Industrial Park, 518057 Shenzhen, China Made in China IVK 

Other labels

Patient tables	 SIEMENS TRUMPF Medical Systeme GmbH • Co. KG Carl-Zeiss-Strasse 1-3 • 73439 Göttingen Made in Germany TIM Dockable Table D21 3T REF: 4756354 SN: 12811001 Net Wt: 240kg 2012-10 IVK TP1 MODEL No. 10849723 (S) Serial No. 12811001
----------------	--

	
Charging station	
Respiratory cushion with pressure hose	 
Respiratory belt	

	
Application cushion (Cushion PERU 098)	
Instrument table, non-magnetic.	
Physiological Measurement Unit (PMU)	Pulse unit:  ECG and respiratory unit:

Description of the symbols on the labels

GTIN	Global GTIN Item Number
Model	Model
	Serial number
	Manufacturer
	Catalog Number
	Do not throw away. Take it to a recycling point.
	Product compliance with European regulations
	Authorized Representative in the European Community
	Caution! Refer to the instructions for use.
	Refer to the instructions for use.

	Class item II
	Working part type BF
	Working part type B

6 Sterility

MR system device is not intended to be sterile.

7 The list of International Standards and Directives applied to the Medical Device



Reference No	Title of the standard
IEC 60601-1:2012 EN 60601-1:2006 + Cor.:2012 + A1:2013	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
IEC 60601-2-33: 2010/A1:2013 /A2:2015 EN 60601-2- 33:2010/A2:2015	Medical electrical equipment - Part 2-33: Particular requirements for the basic safety and essential performance of magnetic resonance equipment for medical diagnosis
IEC 60601-1-2:2014 EN 60601-1-2:2015	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances - Requirements and tests
IEC 60601-1- 6:2010+A1:2013 EN 60601-1- 6:2010+A1:2015	Medical electrical equipment – Part 1-6: General requirements for safety – Collateral standard: Usability
IEC 60825-1:2014 EN 60825-1:2014	Safety of laser products - Part 1: Equipment classification, requirements and user's guide
IEC 62304:2006+A1:2015 EN 62304:2006 + Cor.: 2008+ A1:2015	Medical device software - Software lifecycle processes
IEC 62366:2015 EN 62366:2015	Medical devices - Application of usability engineering to medical devices
EN 1041:2008/A1:2013	Information supplied by the manufacturer of medical devices
EN ISO 15223-1:2016	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements

ISO 10993-1:2018-08 EN ISO 10993-1:2020-12	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
ISO 14971:2007 EN ISO 14971:2012	Medical devices - Application of risk management to medical devices
ISO 13485:2003+AC:2007 EN ISO 13485:2012	Medical devices - Quality management systems - Requirements for regulatory purposes

8 List of equipment and measurement tools necessary for quality control of the medical device











The following accessories for quality assurance are available at the customer. The phantoms are not calibrated.

	Plastic bottle, 15ml Solutes per 1000 g H2O dist.: 3g MnCl2 * 4 H2O + 5g NaCl
	Plastic bottle, 1900ml Solutes per 1000 g H2O dist.: 3.75g NiSO4* 6 H2O + 5g NaCl
	Plastic bottle, 5300ml Solutes per 1000 g H2O dist.: 3.75g NiSO4* 6 H2O + 5g NaCl
	Spherical phantom D165 Diameter: 165 mm Contents: 2570 ml Solutes per 1000 g H2O dist.: 1.25g NiSO4 * 6 H2O

		Spherical phantom D240 Diameter: 240 mm Contents: 7300 ml Solutes: MARCOL-Oil + 0.011g MACROLEX blue
		Hand/Wrist 16 phantom Contents: 1100 ml Solutes per 1000 g H2O dist.: 3.75g NiSO4* 6 H2O + 5g NaCl
		Foot/Ankle 16 phantom Contents: 3300 ml Solutes per 1000 g H2O dist.: 3.75g NiSO4* 6 H2O + 5g NaCl
		Quality assurance coil for 3T Endorectal
		Special phantom for Endorectal A powder is included with the special phantom for the Endorectal coil. It is used to prepare the phantom fluid.

Phantoms Glossary – see Appendix 1

The following phantom holders need to be used for the Local Coil's quality assurance

		
Phantom holder 17	Phantom holder 18	Phantom holder 23
		
Phantom holder 24	Phantom holder 25	Phantom holder 26
		
Phantom holder 27	Phantom holder 28	Phantom holder 29
		
Phantom holder 31	Phantom holder 32	Phantom holder 40
		
Phantom holder 45	Endorectal coil simulator	Holder for spheric phantom

 Phantom holder Breast Sen		
--	--	--

Phantom holders Glossary – see Appendix 1

9 Transportation requirements and conditions

Transportation requirements

- The magnet can be shipped and transported worldwide with or without liquid helium. Maximum possible helium level during shipping: 95%.
- The magnet is classified as “Dangerous goods” as defined by the ADR, IMDG and IATA codes and their requirements. Liquid helium is classified under UN1963 (Helium, refrigerated liquid).
- Transportation of the system containing liquid helium must be in accordance with the appropriate national and international codes.
- Using an air-ride truck with or without a shipping pallet the following maximum peak acceleration measured at the magnet mounting points during transportation shall not exceed:
 - 5.0g vertical shock load applied to the cryostat mounting points.
 - 2.0g horizontal shock load applied to the cryostat mounting points in any direction.
 - Shock is defined as a dynamic force with duration of no more than 0.5 seconds.
- Road transport is allowed up to 15.000 miles.
- Road transport may be with or without an approved shipping pallet. Without NO shipping pallet, the magnet must be supported using transport beams.
- Road vehicles must be fitted with air ride suspension, whether or not a shipping pallet is used.
- Always transport a magnet using vehicles fitted with air ride suspension.
- Road vehicles must be suitable for the intended usage, be properly maintained at all times, and comply with all local laws and regulations.
- Relevant requirements of the ADR code shall be satisfied.
- The vehicle and operator must be approved by Siemens.
- This specification does not constitute acceptance for use as part of a mobile system. Mobile magnets require a specific trailer design for continuous road transportation.
- Air transport is allowed up to 30.000 miles.

- If air transport is used for any section of the route, an approved shipping pallet must be used for the entire transportation process.
- Relevant requirements of the IATA code must be satisfied.
- Sea transport is allowed up to 30.000 miles.
- If sea shipment is used for any section of the route, an approved (magnet system related) shipping pallet must be used for the entire transportation process.
- Relevant requirements of the IMDG code must be satisfied.
- Rail Transportation is not allowed!
- Transportation through the channel tunnel (UK to France) is permissible if the system is mounted on an approved air ride truck. The air ride suspension may be inactive.
- Do not use a fork lift for lifting and transporting the magnet!
- Lifting the magnet must always be conducted by personnel who are trained in handling components and medical equipment of this size and weight.
- Lifting during transportation including the final movement onto the installed site must be performed using a suitable sized crane capable of lifting the full magnet weight.
- Suitable straps or chains capable of lifting the full magnet weight must be used. The straps or chains must be secured to the magnet lifting points using the four shackles supplied and fitted to the magnet.
- Use four suitable rated slings or chains (minimum of 3 m in length) for lifting the full system weight.
- Check the magnet components are not fouled by the lifting straps or chains.
- Do not stand on the Magnet. Damage to components may occur.
- All transport activities including delivery of the magnet and the system components to the MR suite may be performed by a licensed carrier only.

Transportation specifications – Magnet

Transport specifications for the Magnet	
Ambient temperature	-20°C to +60°C
Relative humidity	10% to 95% non-condensing (must be kept dry if unpacked)
Atmospheric pressure	20 kPa to 106 kPa

Transportation specifications – System Components

Transport specifications for the system components	
Ambient temperature	-25°C to +60°C
Relative humidity	10% to 95%
Atmospheric pressure	70 kPa to 105 kPa

10 Shelf life

Siemens considers the subject devices exempt from having a shelf life due to the fact that the devices are manufactured upon order, and are therefore not susceptible to degradation that would lead to functional failure.

Once an order is placed and a machine is manufactured, it is transported, installed and QC and tune up tests are completed on site to ensure proper function of the MR scanner. The MR scanners do not spend time on the "shelf" that would pose a risk to the performance of the device or its components.

The expected service life of 10 years of operation under the environmental conditions indicated in the instructions for use as well as with all preventive and corrective maintenance performed by qualified technical personnel authorized by Siemens Healthcare GmbH.

Refurbished systems are considered not to be different compared to standard (non-refurbished) systems regarding the expected service life, since production- and test-procedures including service and maintenance are identical to systems which have not been refurbished.

If the systems are used beyond this time, additional maintenance and service may be required in order to maintain system performance and reliability.

Shelf life of ECG electrodes, non-magnetic, disposable: 24 months from date of manufacture

10.1 Storage specifications – Magnet

Storage specifications for the Magnet	
Ambient temperature	-20°C to +60°C
Relative humidity	10% to 90% non-condensing (must be kept dry if unpacked)
Atmospheric pressure	250 mbar to 1060 mbar

10.2 Storage specifications – System Components

Storage specifications for the system components	
Ambient temperature	-20°C to +60°C
Relative humidity	10% to 90%
Atmospheric pressure	250 mbar to 1060 mbar

11 Storage and Shelf Life

Siemens considers the subject devices exempt from having a shelf life due to the fact that the devices are manufactured upon order, and are therefore not susceptible to degradation that would lead to functional failure.

Once an order is placed and a machine is manufactured, it is transported, installed and QC and tune up tests are completed on site to ensure proper function of the MR scanner. The MR scanners do not spend time on the "shelf" that would pose a risk to the performance of the device or its components.

The expected service life of 10 years of operation under the environmental conditions indicated in the instructions for use as well as with all preventive and corrective maintenance performed by qualified technical personnel authorized by Siemens Healthcare GmbH.

Refurbished systems are considered not to be different compared to standard (non-refurbished) systems regarding the expected service life, since production- and test-procedures including service and maintenance are identical to systems which have not been refurbished.

If the systems are used beyond this time, additional maintenance and service may be required in order to maintain system performance and reliability.

12 Usage conditions

12.1 Installation requirements

For the following topics exist installation procedures, described in detail in the installation instruction.

General information

- Introduction
- Requirements
- Overview of the Components / Items included in the Shipment
- Procedural flowchart

Examination Room

- General Information
- Patient Table Base Plate
- Magnet
- Patient Table
- Magnet Stop button

Equipment Room

- General Information
- GPA/EPC
- RF-filter plate
- SEP with Helium Compressor (optional)
- Helium Compressor

Operator's Room

- General Information
- syngo Acq. WP (host computer)
- Console Components
- Installing the alarm box
- Wireless PMU Sensor Kit
- MR Workplace table (optional)
- Remote Monitor (optional)
- syngo.via

Establishing Water and Helium Connections

- General Information
- Water connections
- Helium Pressure Lines

Establishing System Connections

- General Information
- Wiring Diagrams
- On-Site Voltage Connection
- Protective Earth Connection
- Protective Earth Connection Test
- Line Voltage Cable
- Gradient Cables
- Ramp cables
- RF Connection
- Patient Table Connection
- Fiber-Optic Cable
- Door Switch
- fMRI Trigger Converter
- Remote Alarm Contact
- Gradient Supervision
- Footswitch (optional)
- Comfort Kit (optional)
- Patient Video (optional)

Line Voltage/System Cooling

- General Information
- Line Voltage Adaptation
- System Cooling
- Test the Magnet Stop buttons (ERDU)
- Status LEDs on the Alarm Box (Troubleshooting)

Final Work Steps

- General Information
- Removing the Eye Bolts and the Transport Device
- Positioning the Container for dripping Water
- System Covers
- Cleaning the System
- Returning Installation Aids
- Report
- Waste Disposal

12.2 Information about maintenance (general and on the territory of Russia)

The system owner manual describes the maintenance activities which need to be performed by the system owner or by service personnel.

With Siemens Service or service personnel, a group of specially trained persons who are authorized by Siemens to perform certain maintenance activities is meant. References to "Siemens Service" include service personnel authorized by Siemens.

Maintenance at regular intervals

In the interest of the safety of patients, operating personnel, and third parties, it is strongly recommended that only authorized personnel perform the maintenance procedures provided by Siemens. System checks should be conducted more frequently if the system is operated under extreme conditions. Within the system owner manuals also a maintenance plan is provided.

Daily checks

In the course of system operation, technical and constructional changes may have been made to the MR system and its environment. It must be ensured that all components function satisfactorily and that the changes do not present hazardous conditions.

After system acceptance by the customer, a daily visual inspection should be performed for the MR system with respect to the following constructional changes:

- Changes in the environment at the output of the exhaust line (e.g. window installed retrospectively, inlet and outlet of air-conditioning units, new buildings or temporarily installed containers)
- Changes to the air conditioning unit or venting system (e.g. by adding air inlets and outlets in neighboring rooms)
- Installation of additional MR systems (e.g. prohibited use of the same exhaust line for several MR systems)
- Constructional changes inside and outside the examination room

The examination room must be checked for (newly introduced) magnetic parts.

Annual checks

The annual technical safety inspections are listed in this system owner manual and may only be performed by Siemens Service.

Siemens Healthcare GmbH provides this documentation "as is" without the assumption of any liability under any theory of law.

The content described herein requires superior understanding of our equipment and may only be performed by qualified personnel who are specially trained for such installation and/or service.

12.3 MAGNETOM Maintenance plan

The maintenance plan shows in tabular form the activities that have to be performed by qualified and authorized service engineers in the course of periodic maintenance work.

The maintenance work is subdivided into:

- Safety check (including safety-related tests)
- Preventive maintenance
- Quality and function tests

Detailed working instructions for all maintenance work are provided in the service documentation for MAGNETOM Prisma.

The RF-Cabin and the door of the RF-Cabin are not part of the MAGNETOM system maintenance and not considered in these instructions. The RF-Cabin and the door of the RF-Cabin door have to be maintained according to the instructions of the manufacturer.

Safety check (including safety-related tests)

The following checks are intended to ensure the safety of the system. Where appropriate, preventive measures have to be adopted or repairs performed. The points to be checked are generally regulated by laws and standards. The Safety-related tests according to §6 MPBetreibV (Germany) are mandatory. The specified checks and intervals correspond to the minimum requirements. Compliance with stricter national legislation may occasionally be necessary.

Object or function	Reason	What is checked	Interval
General Tests			
System	Safety of patient, personnel, and system	Visual inspection of the system	Annually
Visible cabling, cable routing	Safety of patient, personnel, and system	Visual inspection of cables, cable routing	Annually
Phantom	Safety of patient, personnel, system, and environment	Phantoms cannot have any defects or air bubbles. Replace or refill as necessary.	Annually
Options	Safety of patient, personnel, and System	Visual inspection of options	Annually
Accessories	Safety of patient, personnel, and system	Visual inspection of accessories (According to manufacturer recommendations)	Annually
Magnet	Optimum function within specifications Safety of patient, personnel, and system	Check pressure of the magnet Check for ice formation Check for leaks	Annually
Patient fan	Safety of patient and personnel	Check / replace the filter	Annually
Quench system	Safety of patient, personnel, and system	Visual inspection of the quench tube	Annually

Object or function	Reason	What is checked	Interval
		- Check for water in the quench valve - Check quench tube and isolation - Check quench tube outlet for obstructions - Check outlet rain cover - Check protective mesh - Check quench tube outlet, restricted area and warning signs	
Operating manuals	Safety of patient, personnel, and system	Required operating manuals are available and legible	Annually
User icons, button labeling, warning labels, warning signs	Safety of patient, personnel, and system	Visual inspection of the user icons, button labeling, and warning labels	Annually
Door to the examination room	Safety of patient and personnel	Check door lock	Annually
Exclusion zone for magnetic field	Prevent hazards caused by magnetic field	The 0.5 mT zone is marked The examination room is identified with the warning "strong magnetic field"	Annually
Electrical tests			
Protective conductor MRSC (MRWP)	Protecting personnel from electric shocks. (The protective conductor resistance can change during system operation, i.e. through oxidation, corrosion, and cable breaks.)	Protective conductor resistance of the MRSC (MRWP)	Annually
Protective conductor Tim Dockable Table	Protecting personnel from electric shocks. (The protective conductor resistance can change during system operation, i.e. through oxidation, corrosion, and cable breaks.)	Protective conductor resistance of the Tim Dockable Table	Annually
Protective conductor system	Protecting patients and personnel from electrical shock. (The protective conductor resistance can change during system	Protective conductor resistance of the entire system	Every 2 years

Object or function	Reason	What is checked	Interval
	operation, i.e. through oxidation, corrosion, and cable breaks.)		
Function			
Emergency power off circuit (if installed)	Safety of patient, personnel, and system	Checking function of the power circuit	Annually
Fix patient table (option)	Safety of patient, personnel, and system	Inspect table movement and end switches Operation of the emergency stop and safety switches Operation of the emergency release • Checking tolerance limits • Marking of the hazard area	Annually
Tim Dockable Table (option)	Safety of patient and personnel	Inspect table movement and end switches Operation of the emergency stop and safety switches Operation of the emergency release • Checking tolerance limits • Marking of the hazard area • Operation of the locking mechanism • Inspect the emergency undocking function • Checking the chassis and the wheels	Annually
Intercom, pneumatic (squeeze) bulb	Communication with the patient	Function of intercom Function of pneumatic (squeeze) bulb	Annually
Monitoring the magnet (MSUP)	Prevent hazards caused by magnet	Operation of the magnet supervision and the magnet stop	Annually
Magnet stop function (ERDU)			
Gradient supervision	Optimum function within specifications	Function of the gradient supervision	Annually

Object or function	Reason	What is checked	Interval
	Safety of patient, personnel, and system		
QA (Quality Assurance) measurements	Optimum function within specifications Patient protection	Interaction of various system components according to the guaranteed system characteristics using the complete QA measurement	Annually

Leakage Current Measurements

National regulations have to be observed, e.g. IEC 62353/DIN EN 62353 in Germany:

- Instructions for Measuring Equipment Leakage Current
According to IEC 62353, DIN EN 62353, it is not necessary to measure the equipment leakage current. Reasons:
 - the protective conductor is permanently connected
 - compliance with the on-site protective standards, e.g. DIN VDE 0100-710, are described in the planning documents. The system owner has to comply with the on-site protective standards.
 - The protective conductor is checked at regular intervals and if required after repair and service-related activities.
- Instructions for Measuring the Patient Leakage Current at the Patient Table (type B part is applied)
It is not necessary to measure the patient leakage current at the patient table. Reasons:
 - applied part has a non-conducting surface,
 - live parts are separated by metallic components with a protective conductor connection, as per IEC 60601-1, section 17.a.2 (2nd edition), and section 8.5.2 (3rd edition).
- Instructions for measuring the Patient Leakage Current at the MR Local Coils (type B or BF part is applied)
It is not necessary to measure the patient leakage current at the local coils. Reasons:
 - applied part has a non-conducting surface,
 - the live parts are separated from the applied part by a power supply with double insulation (2 MOPP),
 - applied part does not have dangerous voltages, as per IEC 60601-1, section 17.a.3 and 17.a.4 (2nd edition), and section 8.5.2 (3rd edition).
- Instructions for measuring the Patient Leakage Current at the Body Coil (type B part applied)
It is not necessary to measure the patient leakage current at the body coil. Reasons:
 - applied part has a non-conducting surface,

- the live parts are separated from the applied part by a power supply/transformer with d Double insulation (2 MOPP) as per IEC 60601-1, section 17.a.3 and 17.a.4 (2nd edition), and section 8.5.2 (3rd edition).
- Instructions for measuring the Patient Leakage Current at the ECG/Pulse Module (type BF part applied)
It is not necessary to measure the patient leakage current at the applied part of the ECG/Pulse module. Reasons:
 - the applied part of the ECG/Pulse module is powered by a battery,
 - the applied part of the ECG/Pulse module is completely electrically isolated by a wireless RF operation.

Preventive maintenance

The purpose of preventive maintenance is to keep unforeseen failures to a minimum. This satisfies the prerequisites for system compliance with the guaranteed, long-term characteristics. The effects of different operating conditions (full or partial load operation, temperature, size of dust particles, humidity, gases, vapors) are checked and the condition of parts subject to wear is determined by recording and analyzing characteristic values. Preventive measures must be adopted or repairs must be undertaken as appropriate. The specified maintenance intervals correspond to the minimum requirements. Compliance with stricter national legislation may occasionally be necessary.

Object or function	Reason	What is checked	Interval
Cooling System	Preventive measure to avoid overheating	• Check for leaks and condensation • Check/replace water filter • Check water pressure and add water if necessary	Annually*
Air filter, fans	Preventive measure to avoid pollution	• Replace air filter • Check functionality of fans	Annually*
Chiller (option)	Preventive measure to avoid overheating	According to manufacturer recommendations	
Color LCD monitor (1)	Optimum function within specifications	Check/adjust according to specifications	Annually*
Examination room door	Preventive measure to avoid wear and tear	According to manufacturer recommendations	
Magnet (OR65)	Optimum function within specifications	Check/set pressure on magnet and compressor	Annually
	Safety of patient, personnel, and system	LHe Refill	n.a. (2)
		Replace the cold head	n.a. (3)
Helium compressor	Optimum function within specifications	Check for leaks	Annually
Adsorber	Optimum function within specifications	Replacing the adsorber	Every 36 months

Object or function	Reason	What is checked	Interval
Fix Patient Table	Preventive measure to avoid wear and tear	• Check the hydraulic system • Check/refill hydraulic oil	Annually
Tim Dockable Table (option)	Preventive measure to avoid wear and tear	• Check the docking mechanism • Check the hydraulic system • Check/refill hydraulic oil • Check/clean connectors of the docking station •	Annually*
Comfort kit (option)	Preventive measure to avoid overheating	Replace air filter	Every 24 months
Software	Preventive measure against data loss in the case of hard disk failure	• Save dynamic data • Clean up directories	Annually

* For these components the initial maintenance interval is 18 months. The recurrent maintenance interval is every 12 months (annually).

- (1) A standard monitor without calibration is not suitable for diagnostic purposes. Please consider the initial acceptance testing for image display devices and the follow-up service for constancy testing on a regular base, as offered by Siemens service.
- (2) Zero helium boil off rate for typical clinical use, depending on sequence and operating time with running helium compressor. The system needs to be serviced at regular intervals. Undisturbed magnet cooling for 24 hours and 7 days a week.
- (3) If the system is not connected to Siemens Remote Service (SRS) the replacement has to be performed according to manufacturer recommendations. The replacement is performed "on demand" when the system is connected to Siemens Remote Service (SRS).

Quality and function tests

Quality and function tests are used to check whether the system complies with the guaranteed characteristics. Image quality tests determine deviations from the original status. In the case of deviations, preventive measures have to be adopted or repairs performed as appropriate.

The specified maintenance intervals correspond to the minimum requirements. Compliance with stricter national legislation may occasionally be necessary.

Object or function	Reason	What is checked	Interval
QA measurements	Optimum function within specifications	Interaction of all system components according to the guaranteed characteristics	QA has to be performed during safety related tests

Service Contacts are:

	Siemens Healthcare GmbH – Customer Service (CS)	Siemens Healthcare LLC
Address:	Henkestr. 127, 91052 Erlangen, Germany	Dubininskaya Str. 96, 115093 Moscow, Russia
Telephone:	+49 9131 84-0	+7 495 737 12 52
Web:	www.siemens-healthineers.com	www.siemens-healthineers.com/ru

12.4 Information about calibration of the MD / requirements to calibration (if applicable)

There is no calibration required for the MR systems. Nevertheless, there are different checks and tests which have to be done as described in the system owner manual.

To ensure correct installation, complete all preliminary steps specified.

For detailed information on preparing the room, transporting the system, and other installation requirements, consult the project manager.

The on-site project manager is responsible for all room and system-related preparations.

The following prerequisites must be met:

- The room must be fully prepared and cleaned.
- Power and water connections must be available.
- The internal resistance of the on-site power supply unit must correspond to the requirements.
- The transport route for the system must be defined and free of obstructions (regarding dimensions, access, load).
- Ensure that the magnet room is properly equipped.
- Ensure that the openings of the RF shield are sufficiently large.
- Ensure that there is sufficient space to install the magnet.
- Take account of narrow passages due to suspended ceilings, carriers, and columns.
- Determine the final location of the magnet in the magnet room.
- Ensure that the primary cold water supply is ready for operation.
- Block off the unloading zone to ensure that passersby are not harmed.
- Do not block escape routes or essential access routes to the hospital.
- The rooms must be properly equipped. The RF-shield, including the opening for inserting the magnet, must be installed. The magnet room must also be equipped with sufficient lighting. The floor must be clear so that the magnet and system components can be transported.
- The four corners of the magnet must be drawn on the floor of the magnet room to determine the magnet's location. Refer to the project plan.
- The floor from the unloading zone to the final location must be covered with metal plates (at least 5 mm thick).

- The tools required for installation must be available on site (steel rollers, hydraulic lifts, etc.).

12.5 Service contacts

Regular inspections should be a part of the annual maintenance performed by Siemens Service under a service contract. We therefore recommend that you take out a service contract. It is also possible to take out a service contract based on system usage. Please contact your authorized local representative in this matter.

The name and address of the local representative in Russia:

Siemens Healthcare LLC
Dubininskaja Str. 96,
115093 Moscow, Russia
Tel.: +7 495 737 12 52
www.siemens-healthineers.com/ru/

Customer's support center (24 hours):

Phone: 8 800 550 50 22
Email: service.ru@siemens-healthineers.com

13 Warranty

The warranty period supposed to be no less than 12 months from the date of commissioning but defined by the contract terms for every purchase.

The system owner is responsible for the software licenses including software warranty conditions.

The specific warranty conditions are usually defined in the contractual agreement between Siemens Healthcare GmbH and the system owner.

Siemens considers the subject devices exempt from having a shelf life due to the fact that the devices are manufactured upon order, and are therefore not susceptible to degradation that would lead to functional failure.

Once an order is placed and a machine is manufactured, it is transported, installed and QC and tune up tests are completed on site to ensure proper function of the MR scanner. The MR scanners do not spend time on the "shelf" that would pose a risk to the performance of the device or its components.

The expected service life of 10 years of operation under the environmental conditions indicated in the instructions for use as well as with all preventive and corrective maintenance performed by qualified technical personnel authorized by Siemens Healthcare GmbH.

Refurbished systems are considered not to be different compared to standard (non-refurbished) systems regarding the expected service life, since production- and test-procedures including service and maintenance are identical to systems which have not been refurbished.

If the systems are used beyond this time, additional maintenance and service may be required in order to maintain system performance and reliability.

14 Contact for reclamation

In case of any questions or reclamations, the following shall be contacted:

Germany	Russia
Siemens Healthcare GmbH – Magnetic Resonance (MR) Henkestr. 127, 91052 Erlangen, Germany Tel.: +49 9131 84-0 www.siemens-healthineers.com	Siemens Healthcare LLC Dubinskaya Str. 96, 115093 Moscow, Russia Tel.: +7 495 737 12 52 Fax: +7 495 737 13 20 www.siemens-healthineers.com/ru/

Phantoms Glossary

The «Phantoms» are intended for calibration, distortion diagnostics and image quality control and are also used when installing various RF-coils. All variants of Phantoms perform the same function, regardless of the name.

Phantoms (no more than 15 pcs., if necessary)	Plastic bottle, 15ml
	Plastic bottle, 1900ml
	Plastic bottle, 5300ml
	Spherical phantom D165
	Spherical phantom D240
	Hand/Wrist 16 phantom
	Foot/Ankle 16 phantom
	Quality assurance coil for 3T Endorectal
	Special phantom for Endorectal

Phantom holders Glossary

All variants of «*Phantom holders*» perform the same function, as they are intended to fix the correct position of the phantoms.

Phantoms holders (no more than 15 pcs., if necessary)	Phantom holder 17
	Phantom holder 18
	Phantom holder 23
	Phantom holder 24
	Phantom holder 25
	Phantom holder 66
	Phantom holder 14
	Phantom holder 27
	Phantom holder 28
	Phantom holder 29
	Phantom holder 31
	Phantom holder 32
	Phantom holder 40
	Phantom holder 45
	Endorectal coil simulator
	Holder for spheric phantom
	Phantom holder Breast Sen

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Herewith, I certify the foregoing signatures,
acknowledged before me, of

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Mrs. Hertle was identified by
presenting her official photo ID.

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TxRx Knee 15

TxRx Knee 15 Flare

Quality Electrodynamics, LLC
6655 Beta Drive, Suite 100
Mayfield Village, OH 44143 USA
Distributed by Siemens Healthcare GmbH

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Quality Electrodynamics, LLC
6655 Beta Drive, Suite 100, Mayfield Village, OH 44143 USA
Distributed by Siemens Healthcare GmbH

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- Head Coil 31P1H

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- Head Coil 23Na1H

Manufactured by
RAPID Biomedical GmbH
Kettelerstr. 3-11, 97222 Rimpf
Germany - www.rapidbiomed.de

³¹P / ¹H
Head Coil 3T

Model Type	V-XQ-HQ-030	
 REF	01487	
 SN	027	
Version	V02	
	2019-11	
Nucleus	³¹ P	¹ H
Frequency	49.9 MHz	123 MHz
Channels	TxRx (quadrature)	TxRx (quadrature)
MR System	Siemens MAGNETOM 3T	

ATTENTION: Do not use if there is damage
of enclosure, connections or cord.

Manufactured by
RAPID Biomedical GmbH
Kaisersstr. 3-11, 97222 Rimpfing
Germany - www.rapidbiomed.de

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- Shoulder Small 16;
- и/или
- Shoulder Large 16

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