*To the Authorities and test laboratories
involved in the product registration in
Russia*

www.brainlab.com

Feldkirchen, September 21, 2015

I here certify that the documents attached are true and correct to the best of my knowledge and belief.

They have been defined accordingly to internal and external requirements and are up to date.

Sincerely,

A handwritten signature in black ink, appearing to read "AS", written over a horizontal line.

Alexander Schwiersch
Regulatory Affairs manager

Alexander.schwiersch@brainlab.com
tel: 089-9915680

AIRO MOBILE CT SYSTEM

User Manual

Copyright 2013, Mobius Imaging LLC. All rights reserved.

TABLE OF CONTENTS

1 GENERAL INFORMATION	1
1.1 Chapter Overview	1
1.1.1 Contents	1
1.2 Contact Data	2
1.2.1 Contacting the Company	2
1.3 User Manual	3
1.3.1 System Described	3
1.3.2 Symbols Used in This Manual	4
1.4 AIRO Mobile CT System	5
1.4.1 Overview	5
1.5 Overview of Components	8
1.5.1 Transport Mode	8
1.5.2 Scan Mode	9
1.6 Legal Information	10
1.6.1 Summary	10
1.7 Training and Documentation	11
1.7.1 Training	11
1.7.2 Applications Training	12
1.7.3 Documentation	13
2 SAFETY	14
2.1 Chapter Overview	14
2.1.1 Contents	14
2.2 Responsibilities	15
2.2.1 After-Sale Operating and Safety Practices	15
2.3 In Case of Emergency	16
2.3.1 Overview	16
2.4 CPR Mode	18
2.4.1 Performing CPR on Patient in System	18
2.5 Draping and Sterile Fields	19
2.5.1 Operating Room Use	19
2.6 Data Safety	20

2.6.1 Protecting AIRO from Computer Viruses.....	20
2.7 Radiation Safety	21
2.7.1 General Radiation Safety Information.....	21
2.7.2 Safety Features Against Excessive Radiation	22
2.7.3 Profile of Stray Radiation	23
2.7.4 Protective Barriers.....	25
2.7.5 Protection by Distance.....	26
2.7.6 Assessing Patient Risk.....	27
2.7.7 Reducing Risk By Using Preset Dosages	28
2.7.8 Laser Safety.....	29
2.7.9 Reducing Risk from Incorrect Scan.....	30
2.7.10 Reduce Risk of Skin Burns.....	31
2.8 Mechanical Motion Safety	32
2.8.1 Overview.....	32
2.8.2 Patient Environment	33
2.8.3 Safety Precautions in the Patient Environment	34
2.9 Electrical Safety.....	35
2.9.1 Overview.....	35
2.10 Transportation Safety.....	36
2.10.1 Overview	36
2.11 Fire Safety.....	37
2.11.1 Fire Extinguishers	37
2.12 Power Failure	38
2.12.1 In Case of Power Failure	38
2.12.2 Restoring a Session	40
3 TRANSPORT MODE.....	41
3.1 Chapter Overview	41
3.1.1 Contents	41
3.2 Controls	42
3.2.1 Pendant and Power	42
3.2.2 Power Switches.....	43
3.2.3 Power and Network Cables.....	44
3.3 Transportation.....	45
3.3.1 Turn On Power	45
3.3.2 Starting from Parked Position.....	46
3.3.3 Starting from Scan Mode	48
3.3.4 Rotating and Latching the Ring.....	49
3.3.5 Raising the System for Transport	50
3.3.6 Transportation.....	51
3.3.7 After Transport.....	53

3.3.8 Setup for Scan	54
3.3.9 Lowering the System	56
3.3.10 Unlock and Rotate the Ring	57
3.3.11 Setting the Orientation for Imaging	58
3.3.12 Temporary Storage	59
3.4 Powering Down	60
3.4.1 Set to Standby Mode	60

4 COLUMN AND TABLE

4.1 Chapter Overview	61
4.1.1 Contents	61
4.2 Precautions	62
4.2.1 Overview	62
4.3 Overview of the Column	63
4.3.1 Surgical Column	63
4.3.2 Column	64
4.3.3 Audio Signals	65
4.4 Overview of the Controls	66
4.4.1 Column Keypad and Wireless Remote Control	66
4.4.2 Using the Controls	66
4.5 Column Keypad	67
4.5.1 Overview	67
4.5.2 Column Keypad Visual Signals	70
4.6 Wireless Remote Control	71
4.6.1 Overview	71
4.6.2 Charging the Wireless Remote Control	72
4.6.3 Overview of the Wireless Remote Control	73
4.6.4 Touchscreen	74
4.6.5 Keypad	76
4.7 Controlling the Tabletop	78
4.7.1 Overview	78
4.7.2 Instructions	79
4.7.3 On and Off	80
4.7.4 Lock and Unlock Two-Button Operation on the Column	81
4.7.5 Rotate and Brake the Column	82
4.7.6 Instructions	82
4.7.7 Lift and Lower	83
4.7.8 Zero Position (Transfer Position)	84
4.7.9 Select the Head Position: Column Keypad Only	85
4.7.10 Tilt	86
4.7.11 Trendelenburg and Anti-Trendelenburg	87
4.7.12 Other Tabletop Functions	89

4.8 Functions on Wireless Remote Control Only	90
4.8.1 Functions Specific to Wireless Remote Control	90
4.8.2 Slow	91
4.8.3 Flex and Reflex	92
4.8.4 Light	93
4.8.5 Position	94
4.8.6 Stored Positions	96
4.8.7 Service	97
 5 TRANSFERRING THE TABLE	 100
5.1 Chapter Overview	100
5.1.1 Contents	100
5.2 Transfer Safety	101
5.2.1 Introduction	101
5.2.2 Patient Transfer Safety	102
5.3 Transferring to or from the Shuttle	103
5.3.1 Transferring the Table from the Shuttle to the Column	103
5.3.2 Transferring the Table to the Shuttle	105
5.3.3 Reversing Direction	107
 6 SCANNING	 108
6.1 Chapter Overview	108
6.1.1 Contents	108
6.2 System and Controls	109
6.2.1 AIRO System	109
6.2.2 Pendant Display	110
6.3 Before Scanning	111
6.3.1 Overview	111
6.3.2 Option: Logging In	112
6.3.3 Choosing Scan Mode	113
6.3.4 Connecting System	115
6.4 Moving the Ring in Scan Mode	117
6.4.1 Translating the Ring (z axis)	117
6.5 Prepare Scan Field of View	118
6.5.1 Position and Diameter of Scan Field of View	118
6.5.2 Remove Non-Essential Objects	119
6.5.3 Remove Metal Objects	120
6.6 Welcome Menu	121
6.6.1 Welcome Menu and Functions	121
6.6.2 Menu	123

6.7 Typical Scan Sequence: New Scan	124
6.7.1 New Scan	124
6.7.2 Entering New Patient Data	126
6.7.3 Select Existing Patient Data	127
6.7.4 Entering Patient Data from Hospital Worklist	128
6.7.5 Entering Patient Data from Navigation System	129
6.7.6 Confirming Ring Rotation	130
6.7.7 Selecting Exam Region	131
6.7.8 Selecting Navigated Scan	132
6.7.9 Selecting Patient Orientation and Position	134
6.7.10 Selecting 2D Scout Scan View	135
6.7.11 Setting 2D Scout Scan End Position	136
6.7.12 Setting the 2D Start Position	137
6.7.13 Verifying 2D Parameters and Scanning	138
6.7.14 Selecting 3D Scan Region	139
6.7.15 Move Ring to Start Position	140
6.7.16 Verifying 3D Parameters and Scanning	141
6.7.17 Summary of 3D Image Controls	142
6.7.18 Main Image View	143
6.7.19 Window/Level	144
6.7.20 Image Controls	145
6.7.21 Setting Slice Thickness	146
6.7.22 Sending Scan to DICOM or USB	147
6.7.23 Next Steps	148
6.8 View an Existing Scan	149
6.9 Emergency Scanning	150
6.9.1 Taking Emergency Scans	150
6.10 Editing Patient Data	151
6.10.1 Editing or Adding Patient Data	151
6.11 After Scanning	153
6.11.1 Post-Scanning Steps	153
6.12 Service Functions	154
6.12.1 Service Functions Menu	154
6.12.2 Hospital Service	155
6.12.3 Option: Configure	157
6.12.4 Option: Add New Phantom	158
6.12.5 Option: E-Stop Test	159
6.12.6 Option: Manual Scan	161
6.12.7 Option: Axial Scan	162
6.12.8 Option: Single Shot Scan	163
6.12.9 BrainLab Service	164
7 MAINTENANCE	165
7.1 Chapter Overview	165

7.1.1 Contents	165
7.2 Daily Gain Calibration	166
7.2.1 Overview	166
7.2.2 Preparing for Calibration	167
7.2.3 Perform Gain Calibration	168
7.3 Cleaning, Disinfection, and Sterilization	173
7.3.1 Cleanable Parts	173
7.3.2 Sterilization	175
7.4 Inspection	176
7.4.1 Daily	176
7.5 Preventative Maintenance	177
7.5.1 Preventative Maintenance is Required	177
7.6 Temporary Disuse	178
7.6.1 Temporary Disuse	178
7.7 Disposal Instructions	179
7.7.1 Disposal of System and Batteries	179
 8 QUALITY TESTING	 180
8.1 Chapter Overview	180
8.1.1 Contents	180
8.2 Test Object and Software	181
8.2.1 Description	181
8.3 Testing Image Quality	182
8.3.1 Testing Image Quality with the Phantom	182
8.3.2 Scan a Phantom for Quality Testing	183
8.4 Quality Check Reporting	185
8.4.1 Examining Daily Image Quality Report	185
8.5 Exposure Testing	187
8.5.1 Responsible Organization	187
8.5.2 Measure and Verify CTDI	188
 9 TECHNICAL DESCRIPTION	 190
9.1 Chapter Overview	190
9.1.1 Contents	190

9.2 Introduction	191
9.2.1 Obtaining More Information	191
9.2.2 Warning Regarding Modifications	192
9.2.3 Combining Accessories and Sub-assemblies with AIRO	193
9.3 Technical Description	194
9.3.1 Brief Description and Operating Principle	194
9.3.2 Significant Characteristics	196
9.3.3 How the Device Functions	197
9.3.4 Physical Characteristics	198
9.3.5 Diagrams	199
9.3.6 Pendant Specifications	201
9.3.7 Operating Table Specifications	202
9.3.8 Column Specifications	203
9.3.9 Wireless Remote Control for Column	204
9.3.10 Charger	205
9.4 Special Installation Requirements	206
9.4.1 Overview	206
9.5 Electrical Specifications	207
9.5.1 Specifications	207
9.6 X-ray Specifications	209
9.6.1 Generator (High Voltage Power Supply)	209
9.6.2 X-ray Source Assembly	210
9.6.3 Battery System	211
9.7 Performance Characteristics	212
9.7.1 General Performance Characteristics	212
9.7.2 Range, Accuracy, and Precision	213
9.8 Environmental Requirements	214
9.9 Safety Classifications	216
9.10 Replacing Fuses and Other Parts	217
9.10.1 Replacing Parts	217
9.11 Network/Data Coupling to Other Equipment	218
9.11.1 Overview	218
9.12 Regulatory Compliances and Approvals	219
9.12.1 Summary of Standards Compliance	219
9.12.2 General Electromagnetic Immunity	220
9.12.3 Electromagnetic Immunity	221
9.13 RF Communications Equipment	223
9.13.1 Overview	223
9.14 Symbols	224
9.14.1 Symbols Used	224

9.15 Labels	226
9.15.1 Labels Used	226
9.16 Location of Certification Labels	229
9.16.1 Overview	229
 10 CT USER INFORMATION	 230
10.1 Chapter Overview	230
10.1.1 Contents	230
10.2 No Quick-change X-ray Tubes	231
10.2.1 No Quick-change X-ray Tubes	231
10.3 Dosage Factors and Image Quality	232
10.3.1 Conditions of Operation	232
10.3.2 Dose Information	233
10.3.3 Statement of Maximum Deviations	239
10.3.4 Imaging Performance Information	240
10.3.5 Technical Statements	243
10.3.6 Definitions of the Measurement Bases of the Technique Factors	247
10.3.7 60601 Requirements: CT USERS	248
10.4 Deterministic Effects	249
10.4.1 Radiation Dose Levels	249
10.4.2 Prevailing Radiation Dose	250
 11 ALLOWED EQUIPMENT	 251
11.1 Chapter Overview	251
11.1.1 Contents	251
11.2 Compatibility with Other Medical Devices and Software	252
11.2.1 Medical Devices	252
11.2.2 Medical Software	253
11.2.3 Accessories List	254
 12 TROUBLESHOOTING	 255
12.1 Chapter Overview	255
12.1.1 Contents	255
12.2 Troubleshooting	256
12.2.1 Problems and Suggested Actions	256
12.3 Metal Artifacts in Scans	258
12.3.1 Recognizing Metal Artifacts in CT Scans	258

12.4 Transport Failure	259
12.4.1 Transport Failure	259
12.4.2 Base is Lifted but System Will Not Move	260
12.4.3 Cannot Raise System	261
12.5 Table Transfer Errors	263
12.5.1 Transfer Errors	263
12.6 Table Malfunctions	264
12.6.1 Situation: Unable to Remove Table	264
12.6.2 Column Fails to Reach Zero Position	265
12.6.3 Table Will Not Release	267
12.7 Malfunctions and Return Instructions	268
12.7.1 Malfunctions	268
12.7.2 Return Instructions	269

1 GENERAL INFORMATION

1.1 Chapter Overview

1.1.1 Contents

Topics Covered

Section	See
Contact Data	1.2
User Manual	1.3
AIRO Mobile CT System	1.4
Overview of Components	1.5
Legal Information	1.6
Training and Documentation	1.7

1.2 Contact Data

1.2.1 Contacting the Company

Support

If you cannot find information you need in this manual, or if you have questions or problems, contact Brainlab customer support:

Region	Telephone and Fax	Email
United States and Canada	Tel: (800) 597-5911 Fax: (708) 409-1619	us.support@brainlab.com
Europe	Tel: +49 89 991568-44 Fax: +49 89 991568-811	support@brainlab.com
France and French-speaking regions	Tel: +33 800 676 030	

Distributor and Authorized EU Representative

Brainlab AG
Kapellenstr. 12
85622 Feldkirchen
Germany

www.brainlab.com
Tel: +49 89 99 15 68 0
Fax: +49 89 99 15 68 33
contact@brainlab.com

Manufacturer

Mobius Imaging LLC
323 W. Main Street
Ayer MA 01432
USA

www.mobiusimaging.com
Tel: +1 (978) 796-5068
info@mobiusimaging.com

Authorized representative in Russia

Neuroproject Ltd., 127486, Moscow
Korovinskoe road, 10 building 2, office 43
cell/fax: +7(495)660 -8173
e-mail : info@neurop.ru
web: neurop.ru

1.3 User Manual

1.3.1 System Described

Catalog Number
MobiCT-32

Configuration:

AIRO Mobile CT System with accessories

I. AIRO Mobile CT System:

2. Base
3. Gimbal
4. Power panel
5. Ring
6. Pendant Display
7. Casters – 4pcs
8. Front View Camera
9. Gammex phantom for calibration with case and platform
10. Software for automatic transfer and registration AIRO images for cranial and spine cases

II. Accessories:

1. Batteries
2. X-ray tube

1.3.2 Symbols Used in This Manual

General	Warnings and cautions appear in this document and on the equipment. Heed all warnings and cautions when using AIRO .

Warnings

Warnings are indicated by triangular warning symbols. They contain safety-critical information regarding possible injury, death or other serious consequences associated with equipment misuse.

Cautions

Cautions are indicated by circular caution symbols. They contain safety-critical information regarding possible problems with the device. Such problems include device malfunctions, device failure, damage to device or damage to property.

Notes	<i>NOTE: Notes are formatted in italic type and indicate additional useful hints.</i>

Commands	The names of keys ("buttons") and icons are in bold.

1.4 AIRO Mobile CT System

1.4.1 Overview

Device Description **AIRO Mobile CT System** ("AIRO") is a mobile computed tomography (CT) system. AIRO includes a mobile base that supports the gimbal, ring, a surgical column, and controls.

See 9 *TECHNICAL DESCRIPTION* for more detailed information.


Figure 1

**Catalog
Number**

MobiCT-32

Indications for Use AIRO® is intended for use in several environments such as mobile or fixed General Radiology, Intensive Care Unit (ICU), Emergency Department (ED), Surgical/Operating Room (OR), Clinic or Office, to be used for non-contrast X-ray computed tomography applications for anatomy that can be imaged in the 107cm aperture excluding pediatric patients and patients weighing over 400 lbs (182 kg). Airo CT System is intended for use in several environments such as Mobile or fixed General Radiology, Intensive Care Unit (ICU), Emergency Department (ED), Surgical/Operating Room (OR), Clinic or Office. In addition, both Systems provide non-contrast CT images that assist in a range of clinical applications which includes Neuro (cranial and spinal), ENT, head and neck, orthopedic and general surgery.

Contraindications AIRO® is contraindicated as the principal means of guidance during surgical procedures. The physician must verify navigation accuracy using known anatomical landmarks or an equivalent verification method when used in surgical procedure.

Intended Use AIRO is intended for use in several environments such as mobile or fixed General Radiology, Intensive Care Unit (ICU), Emergency Department (ED), Surgical/Operating Room (OR), Clinic or Office. AIRO will provide non-contrast CT images that assist in a range of clinical applications which includes neuro (cranial and spinal), ENT, head and neck, orthopedic and general surgery.

This device is not for use on pediatric patients (defined as birth to 21 years of age) which includes patients less than approximately 50 kg (110 lbs) in weight and 150 cm (59 in) in height; these height and weight measurements approximately correspond to that of an average 12 year old or a 5th percentile U.S. adult female.

Medical Purpose AIRO can provide an aid to diagnosis when used by a qualified physician and is intended for computed tomography applications that can be imaged in the 107cm aperture.

AIRO is not designed to be relied on as the principal means of guidance during invasive procedures when real time imaging is needed. (e.g., CT Fluoro). Risk of patient injury.

The physician must verify navigation accuracy using known anatomical landmarks or an equivalent verification method.

Side effect There are no Side effect after scanning in AIRO.

Proposed clinical

use

Diagnosis of body parts:

- brain, its structure, eyes, inner ear, sinuses;
- neck, shoulders, cervical spine and blood vessels of the neck;
- chest cavity, heart, aorta, lungs, mediastinum;
- upper stomach, liver, kidney, pancreas and other organs of the abdominal cavity;
- pelvis, hip, male and female reproductive system, bladder;

AIRO system is use for the control and conduct of medical and diagnostic procedures in the operative field.

Intended Users and Operators AIRO is intended to be operated by medical personnel trained in its operation by authorized trainers. CT functions may be further limited by local regulations.

AIRO is only intended to be shipped or serviced by qualified Brainlab service personnel and manufacturing technicians.

Only trained medical personnel may operate system components and accessory instrumentation.

See 1.7 *Training and Documentation*.

Patient Group This device is to be used at the request of the supervising physician on patients weighing less than 182 kg (400 lbs) and over 21 years of age.

This instrument is designed to image patients with a mass of less than 182 kg (400 lbs). For patients heavier than 182 kg, device X-ray exposure limitations result in poor image quality and unused dose.

Refer to the instruction manual for the operating table you will be using to determine the respective effective load limits.

Do not exceed the indicated load limit of the operating tabletop. Failure to adhere to proper loading limits can lead to patient injury due to tabletop structural integrity failure.

Intended Location **AIRO** is to be used only in areas that the responsible organization has designated and declared suitable for safe operation.

Place of Control AIRO is controlled inside the treatment room during scan mode. For scanning purposes, **AIRO** is not intended to be controlled from outside the treatment room. The operator of the equipment needs to have a direct and clear line of sight to the device and to the patient when in use and when moving the equipment.

Tolerance to Water The moisture tolerance of the device is component dependent. Nothing more wet than a wrung out drip-less towel should be used anywhere near the device. The column is rated IPX4. Overall system is IPX0.

Improper Use

Improper use is defined as using AIRO or any part of AIRO in any way inconsistent with the methods and precautions described in this manual, or with the methods described in authorized Brainlab training, or use by anyone who is not properly trained and medically qualified.

1.5 Overview of Components

1.5.1 Transport Mode

Diagram: Transport Mode

Figure 2

No.	Component	Function/Usage
1	Ring	Scans patient and sends data to Pendant.
2	Pendant Display	Controls system movement, imaging, and data storage.
3	Power Panel	Key switch and On button, USB port.
4	Gimbal	Holds the ring and Pendant Display.
5	Base	Supports and moves the system. Power and network connections.
6	Casters	Allow system to be moved. <i>NOTE: Not shown: A mechanically driven main drive wheel located under the base. The main drive wheel powers the motion.</i>
7	Front View Camera	Allows operator to see path on display. Used when moving the system for safe transportation.

1.5.2 Scan Mode

Diagram: Scan
Mode

Figure 3

No.	Component	Function/Usage
1	Ring	Scans patient and sends data to Pendant.
2	Base	Supports and moves the system. Location of power and network connections. <i>NOTE: Casters are retracted into the base.</i>
3	Gimbal	Holds the ring and Pendant Display.
4	Control Panel	Holds key switch, on button, battery indicator lights, and USB port.
5	Pendant Display	Controls system movement, imaging, and data storage.
6	Surgical table Triumph	System AIRO could be connect with surgical tables from TRUMPF Medizin Systeme GmbH +Co. KG

1.6 Legal Information

1.6.1 Summary

Copyright This manual contains proprietary information protected by copyright. No part of this manual may be reproduced or translated without express written permission of Mobius Imaging.

Trademarks • **AIRO®** is a registered trademark of Brainlab AG in the US.

The respective owners reserve all rights to their trademark(s).

Restrictions on Use

US Federal law restricts the sale, distribution and use of this device to, or on the order of, a physician.

This device is limited by US federal law to investigational use for indications not in the “indications for use” statement.

DICOM Conformance Statement A DICOM Conformance Statement is available upon request.

Restriction of Hazardous Substances Directive (RoHS) The AIRO System complies with the Directive on the restriction of the use of certain hazardous substances (RoHS) in electrical and electronic equipment. A RoHS compliance statement is available upon request.

Disposal Instructions

Only dispose of electrical and electronic equipment in accordance with statutory regulations. For information regarding the WEEE (Waste Electrical and Electronic Equipment) directive, visit: <http://www.brainlab.com/weee>. A WEEE Statement is available upon request.

Life Cycle Life Cycle in use is 10 years.

Patents Mobius Imaging, LLC has intellectual property rights relating to technology embodied in the product described in this manual. In particular, and without limitation, these intellectual property rights may include one or more of the patents listed at <http://www.mobiusimaging.com/patent.htm> and one or more additional patents or pending patent applications in the US and other countries.

1.7 Training and Documentation

1.7.1 Training

Required Training Users must receive and complete training from Brainlab or its authorized agents before using AIRO. For training, contact Brainlab.

NOTE: The User Manual does not replace training.

Required Training	Hours
Applications Training	40 hours (5 days)
Refresher Training	As requested

1.7.2 Applications Training

Description Medical personnel that will be operating **AIRO** shall attend a 5 day applications training course consisting of:

- Training Specialist scheduled for training
- Training Presentation (including safety information, basic setup and use of the system)
- User's Manual (hard-copy supplied with system)
- **AIRO** CT System (including all components)
- Hands-on training including transportation and maneuvering the system in a simulated typical environment
- Trumpf Column hands-on training. To include table transfer, operational controls and programming.
- A test at the end of the training course appropriate to the practice/user

NOTE: All other users (non-operators) will be trained in the areas relevant to their interactions with the AIRO CT System. For example, nurses will be trained on patient handling, table controls, and preparation of sterile fields.

Training Materials: Training materials include, but are not limited to:

- User Manual
- AIRO system

1.7.3 Documentation

Intended Audience This User Manual is intended for all members of the clinical team who are trained in the use of the **AIRO** system.

Reading the User Manual This User Manual contains the information necessary to operate **AIRO** in accordance with its specifications and intended use.

It is important that AIRO users:

- Read the user manual carefully as part of training and before attempting to use AIRO
- Always have access to the user manual

Available Manuals

Document	Part Number
AIRO User Manual	MI-42-0001

Other Information Additional information on this product may be available by request. Contact technical support.

№	Parameter description	Unit
1	Overall dimensions in transport position, cm, no more than	60 x 228 x 198
2	Dimensions in the scanning position, cm, no more than	193 x 330 x 190
3	Weight, kg, no more than	975
4	The outer diameter of the gantry, cm, no more than	168
5	The diameter of the aperture, cm, no less than	107
6	The diameter of the image field review, cm, no less than	51,2
7	The width of the gantry, cm, no more than	38
8	Speed during transportation (movement scanner), km / h, no more than	1,6
9	A typical braking distance, cm, no more than	15
10	Moving gantry during scanning cm, not more than	100
11	Turnaround time when scanning, with no more than	1,92
12	The maximum load on the system through the patient support kg, no more than	182
13	Dimensions of the control panel, cm, no more than	26 x 4 x 17
14	Remote control Weight, kg, no more than	1,4
15	Overall dimensions of the touch screen remote control cm, no more than	17 x 17
16	The voltage on the control input, DC, no more than	24
17	Power consumption remote control, W, not more than	5
18	Generator power source and the computer	Battery installation, internal
19	Power supply of battery installation	electricity mains, 1 phase
20	Rated voltage battery power installation, AC	100-240
21	Nominal frequency range power battery installation, Hz	50/60 Гц
22	Fitting the battery power supply, VA, not less than	1500
23	The least favorable power supply, AC, no more	100
24	Current frequency in Hz, no more	50
25	The maximum X-ray tube current at the rated voltage, mA, no more than	250
26	Maximum voltage X-ray tube with a maximum current, kW, not more than	125
27	Voltage X-ray tube current at maximum power, kW	125
28	Voltage X-ray tube current at maximum power mA	250
29	Rated electric power at a voltage of 120 kV, kW, not more than	30
30	Circuit breakers	time-lag
31	Type of fuse	3 AG
32	Current fuse, A, no more than	15
33	Voltage fuses, not more than	250
34	Isolation of power conductors	External (male)
35	The impedance of the power supply circuit, Ohm, not more than	0,6
36	Cable Length cm, not more than	500
37	The cross section of the power cable mm ² , not less than	1,00
38	Source Type of high voltage (generator)	high voltage
39	Location of generator	gantry
40	Generator Cooling Method	Air
41	Rated power kW generator, not less than	32
42	Number of scans using the generator at full power (120 kV × 250 mA × 20 sec.) For an hour at least	2
43	Nominal generator parameters	125 kV, 250 mA, 20 sec. 125 kV, 125 mA, 40 sec.
44	The voltage at the output of the generator, kV range	80-125
45	Voltage deviation,%, not more than	±1
46	Output current, mA range	5–250
47	Deviation current range,%, no more than	±2
48	The maximum average electric power, W, not more than	521
49	Minimal battery power required for normal operation, V, not less than	280

50	The cooling time between the maximum scanning sessions should be, min, not less than	30
51	scanning time, with no more than	42
52	Type X-ray tube	The rotating anode
53	The design of the X-ray tube	The end of the anode grounded to the metal chassis
54	The shape of the beam of X-ray tube	Conical
55	Target X-ray tube	7 "rhenium-tungsten-molybdenum
56	tube voltage kV range	80-125
57	Accuracy of the tube voltage,%, no more than	1
58	Tube current mA, range	5–250
59	The error in the current tube,%, no more than	2
60	The dimensions of the focal spot, mm, not more than	1 x 1
61	Full filtering at maximum voltage (on site and filter beamforming) along the axis (equivalent mm Al at 120 kV), mm, no more than	6,8
62	The heat capacity anode TE mn, at least	2,3
63	The maximum cooling rate of the anode, th. TE / min, not less than	2874
64	Pipe cooling method	Coolant (50% water, 50% propylene glycol) - air
65	Battery battery installation	Li-ion phosphate battery LiFePO4
66	Nominal voltage at the terminals of the battery installation, DC, not less than	358
67	The total energy capacity of the battery installation, watt-hour, no more than	36 080
68	Certificate of characteristics of the battery cell	UL1642 and MEC 62133
69	Certificate of battery installation performance	UN/DOT 38.3
70	Type of scan	spiral
71	The number of slices	32
72	The rotation period, with no more than	1,92
73	Scanning range cm, not more than	100
74	The scanning time of 100 cm with no more	42
75	Pitch, not more	1,415
76	slice thickness, mm, not more than,	1
77	The maximum length of topogram cm*cm, not more than	24
78	Dimensions topogram scanning, cm*cm, not more than	512 x 512
79	Angular position topogram available for scanning, grad	100
80	Dimensions scanning projection radiography cm, not more than	51,2 x 100
81	The angular position available for scanning projection radiography	0° и 90
82	MTF 50%, lp / mm, not less than	0,21
83	MTF 10%, lp / mm, not less than	0,37
84	MTF 2%, lp / mm, not less than	0,6
85	A typical value of the tomographic slice thickness (full width at half maximum), with a length of tomographic 1 mm, mm, not more than	1,5 мм
86	The maximum deflection, mm, not more than	±0,5
87	The physical thickness of the detector, mm, not more than	1.06
88	Full slice collimation, mm, not more than	33,92
89	Connectivity	Gigabit Ethernet, DICOM 3.1
90	Scanning algorithms	Sharp, Standard, Soft
91	A laser positioner, mW, no more than	3,5
92	A laser positioner wavelength nm, no more than	635
93	A laser positioner, safety class	II
94	Typical values CTDVol, head, helical scanning (at 120kV, 240 mA, 326 mAs, slice collimation 1,06mm, width of cut 1.0 mm, feed per revolution 48 mm	64,7

	pitch 1,415, Soft algorithm characteristic time of 200 mm - 10 c) mGy,	
95	Typical values CTDVol, abdomen, helical scanning (at 120kV, 110 mA, 149 mAs, collimation 1,06mm slice, slice width of 1.0 mm feed per revolution 48 mm pitch 1,415, Standard algorithm characteristic time of 200 mm - 10c) mGy	22,3
96	The diameter of the wheel system, cm	15,24
97	The maximum noise level, no more	70 dBA
	Operation conditions	
98	Temperature,	10 – 25 ° C
99	Relative humidity, without condensation %	from 20 to 70
100	Atmospheric pressure, kPa	70 -110
	Conditions of storage and transportation	
101	Temperature ° C	0 - 50
102	Relative humidity, without condensation, %	5-95
103	Atmospheric pressure, kPa	70 -110

2 SAFETY

2.1 Chapter Overview

2.1.1 Contents

Topics Covered

Section	See
Responsibilities	2.2
In Case of Emergency	2.3
CPR Mode	2.4
Draping and Sterile Fields	2.5
Data Safety	2.6
Radiation Safety	2.7
Mechanical Motion Safety	2.8
Electrical Safety	2.9
Transportation Safety	2.10
Fire Safety	2.11
Power Failure	2.12

2.2

2.3 Responsibilities

2.3.1 After-Sale Operating and Safety Practices

Manufacturer The manufacturer assumes no responsibility or liability for after-sale operating and safety practices; nor can it be responsible for personal injury or damage resulting from misuse of its systems.

Owner The owner must designate areas suitable for safe operation and service of the system and ensure that it is used only in these designated areas. It is the responsibility of the owner to ensure that all personnel follow the manufacturer's radiological safety recommendations. The owner is also responsible for ensuring compliance with additional local requirements, regulations, or organizational safety practices such as wearing protective clothing and radiation monitoring devices while using this system, and inspecting equipment.

2.4 In Case of Emergency

2.4.1 Overview

How to Respond to an Emergency


Figure 4

Steps	
1.	Press Emergency Stop a on the System Controls.
2.	Follow emergency safety procedures established by your organization.
3.	Press the Reset Button s .

Emergency Interlock Device

An interface to the Emergency interlock is available upon request.

In order to interface to the Emergency interlock (when installed), the user must provide a normally-closed interlock contact wired between pins seven and eight of an ODU S13L0C-P26MFC0-7200 (or equivalent) connector.

It is recommended that any customer using the Emergency interlock also construct a dongle with the above mentioned connector with the pins shorted directly to one another to enable the system to be transported while not connected to the interlock switch.

**When to Press
Emergency Stop**

Press the emergency stop button:

- If the system moves without a command from the user
- If releasing a key does not stop the movement
- If the ring does not stop as expected
- If the ring moves in the wrong direction during a scan
- If for any reason during scanning or transportation it becomes necessary to stop motion.

If any system malfunctions are detected, immediately shut down the system and contact technical support.

**How to Release the
Patient**

No motorized system movements are possible after an Emergency Stop.

You will be able to push the gimbal and ring out of the way. You will not be able to detach the surgical table from the column.

Steps	
1.	Push the ring away from the pedestal/column as near as possible.
2.	Remove the patient manually, for example, transfer the patient to a gurney.

Malfunction

If you used Emergency Stop because AIRO malfunctioned, do not use the system until it is evaluated by technical support.

2.5 CPR Mode

2.5.1 Performing CPR on Patient in System

Table Position The equipment positions for performing CPR are:


Figure 5

Position of Table for Performing CPR	
a Parallel to base, column in zero position	Acceptable
s Parallel to base, column in zero position	Acceptable
d Any other orientation	Not Acceptable

CPR Options

Option	Action	
Patient is in the bore of the system		Move the ring out of the way. <i>NOTE: You can use the controls or push the gimbal manually to move it out of the way.</i>
The table is not parallel to the base		Rotate the table into a position parallel to the base <i>NOTE: If the table will not rotate, release the table rotation brake. The rotation brake control is on the column keypad and on the wireless remote control.</i>
The table or patient support is not lowered and level		Lower and level the table or patient support

2.6 Draping and Sterile Fields

2.6.1 Operating Room Use

Drape the Patient When AIRO is used in a sterile field environment during intra-operative procedures, add an additional transparent plastic drape over the patient to maintain the sterile environment while scanning. **AIRO** should not be draped.

NOTE: It is important to use transparent (clear) plastic drapes to ensure accurate laser positioning.

2.7 Data Safety

2.7.1 Protecting AIRO from Computer Viruses

Responsibility	The organization using the system is solely responsible for protecting AIRO from computer viruses. <i>See 9.11 Network/Data Coupling to Other Equipment.</i>

2.8 Radiation Safety

2.8.1 General Radiation Safety Information

General Warning

AIRO produces ionizing radiation. Excessive exposure could cause serious injury.

NOTE: Restrict access to AIRO in accordance with local regulations for radiation protection.

Radiation Warning Signals

When X-rays are being produced, AIRO alerts personnel by these means:

- Yellow lights illuminate on the gimbal.
- A yellow light illuminates on the pendant display.
- The pendant display screen says "Scanning."
- AIRO emits a beeping tone.

2.7.2 Safety Features Against Excessive Radiation

How to Stop a Scan


Figure 6

Step
Release Scan a button.

2.7.3 Profile of Stray Radiation

Technique Factors Stray radiation was measured with the body dosimetry phantom (32 cm diameter) placed at isocenter on the patient table, using a Fluke 451P pressurized ion chamber survey meter. Measurements on the horizontal plane were made at isocenter height, and on the vertical plane at the centerline of the table, at 50 cm intervals. Exposures in μR were made at 120 kV, 25 mA, 2 sec (50 mAs), and recalculated to 192 mAs.

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2.7.4 Protective Barriers

How to Minimize Radiation Exposure

The recommended practices for minimizing radiation to staff are:

Personnel	Recommended Location During Scan
CT operator and essential staff	Behind protective barriers (mobile shields) as needed
Other medical support staff	Outside of the room

Specification for Protective Barrier

For personnel remaining in the room, mobile shields, lead aprons, and thyroid collars must be utilized. Mobile barriers should have a leaded glass window. All protective barriers (shields, collars, and aprons) must be rated 0.25 mm Pb equivalent or higher.

The CT operator must be able to view the patient at all times when a scan is performed to be able to act promptly in case of an emergency.

2.7.5 Protection by Distance

Use the Pendant Display

AIRO is designed to be operated from a distance to protect the CT operator from exposure to stray radiation. The CT operator must be able to view the patient at all times when a scan is performed to be able to act promptly in case of a collision.

Step
Use the long cord on the Pendant Display to control the scanner at a distance from behind protective barriers.

2.7.6 Assessing Patient Risk

How to Assess Patient Risk

Steps	See
1. Enter the patient's weight and specify the area to be scanned.	6 SCANNING
2. Review the automatically calculated radiation dosage on the display.	6 SCANNING

NOTE: Always check the displayed projected radiation dose before initiating the scan.

2.7.7 Reducing Risk By Using Preset Dosages

How to Reduce Risk Using Preset Dosages

AIRO offers only one mode of operation: Automatic Exposure Control (AEC). The general aim of an AEC system for CT is to reduce the variation in radiation doses to different patient cross sections. The AIRO uses an AEC system that adjusts the X-ray tube current based upon the size of the patient which is configured by the user through the selection of patient weight, and exam region (body part). A constant tube current is used for an entire examination or scan series.

AIRO uses preset dosages based on the parameters entered for each scan (area, duration, patient weight, etc.).

There is no manual mode with AIRO. As identified in the AIRO Application Guide for Protocols and Principles certain protocols use patient weight as an input to an Automatic Exposure Control (AEC) calculation of constant tube mA. The calculation is as follows:

$$\text{Adapted mA} = \text{Standard mA} * \text{Patient Weight} / \text{Standard Patient Weight}$$

Where:

- Standard mA = mA for body region in a 70 kg patient
- Standard Patient Weight = 70 kg

To reduce the risk of patient overexposure:

Steps	See
1. Enter patient data and scan area.	6 SCANNING
2. Review the automatically calculated radiation dosage on the display.	6 SCANNING
3. Evaluate the result.	

Options	
Dose is appropriate	Initiate scan
Dose is inappropriate	Check patient data and scan selections

NOTE: Using preset dosages does not replace the expertise of the operator.

Implanted and Body-worn Medical Devices

The X-rays produced during scanning may interfere with active implantable medical devices and body worn active medical devices. Contact the manufacturer of the device for information before conducting scans.

Precautions when Defibrillator is Used

There is a potential of patient burns when a defibrillator is used, if proper safety precautions are not observed. These precautions include isolation from:

- Metal Parts
- Conductive pads
- Damp Towels

2.7.8 Laser Safety

How to Use the Positioning Lasers Safely

Be sure to avoid crossing a patient's eyes with the light beam from a positioning laser.

Figure 7

The class II laser modules **a** in the ring emit visible radiation that may be harmful to the human eye.

NOTE: Some users' movement may be restricted during system movement. The laser lights will turn on automatically during scanning to show the progress of the scan beam. Patients may not have the ability to, for example, close their eyes or turn their heads away from the light.

2.7.9 Reducing Risk from Incorrect Scan

How to Reduce Risk From Incorrect Scan

If a scan is initiated with incorrect parameters:


Figure 8

Step
Release Scan a button.

2.7.10 Reduce Risk of Skin Burns

How to Reduce
Risk by Checking
Distance

There is a risk of patient skin burn if patient anatomy is closer than 5 cm (2") to the interior of the system bore during scanning.

Step
Check the source-to-skin distance. Be sure all patient anatomy will remain at least 5 cm (2") away from the inner surface of the ring bore during a scan.

2.8 Mechanical Motion Safety

2.8.1 Overview

How to Avoid Injury from Moving Parts

AIRO moves during normal use. This section covers safety near the moving parts.

NOTE: The person operating the AIRO controls is responsible for evaluating the risk of collisions and avoiding them.

2.8.2 Patient Environment

Definition The Patient Environment is the space in which intentional or unintentional contact can occur between a patient and parts of the system, or between a patient and other persons touching the system.

Dimensions of the Patient Environment (United States and European Union)

Side View:


Figure 9

Top-Down View:


Figure 10

2.8.3 Safety Precautions in the Patient Environment

Safety Precautions

The system is designed to be operated from a distance with the pendant when it is in scan mode. This allows operators to assure patient safety from collisions by observing the motion of the ring and table from multiple viewpoints.

The operator should be able to observe any potential collisions while setting start and end points for the scout scan or with the lasers. The operator can then be confident that the patient will clear the inner bore of the gantry. The operator should observe the patient during the scan (either wearing lead apron or standing behind lead barrier) to abort the scan should the patient move/fall during the scan.

Stand back from moving parts and use the long cord on the pendant display and the wireless remote control to operate the system.

Injuries could result from power driven equipment parts colliding with other moving or stationary items in the environment. To avoid collisions:

- Be sure the patient is placed according to the instructions.
- Clear the surrounding area of equipment that is within range of the motion of the ring.
- Place AIRO so that surgical personnel can move easily around it and access all equipment without collisions.
- Do not place any unnecessary objects on or near AIRO.
- Anticipate the path of movement. Prevent clothing, drapes, IVs, and equipment from getting caught in moving parts.
- Keep continuous visual contact with the system and the patient.

Avoiding Collisions

Situation	Risk Mitigation
Transferring the table	Evaluate the area for possible collisions. Clear the area of objects, furniture, and bystanders before any transfer.
Moving the ring or the operating table	Monitor all movements. Stop movements by releasing the button on the controls.
Positioning tubing, other surgical equipment	Anticipate the system path. Place items so there is no danger of collision or pinching.
Changing table incline	Monitor the changes in patient incline.
Moving the column to zero position	Select Zero Position on the column only after evaluating the area for possible collisions. Pay particular attention to patient incline and any table accessories that are angled downward.

Do Not Lift

There are no parts of AIRO that need to be lifted during authorized use with the exceptions of the pendant display and the remote control.

2.9 Electrical Safety

2.9.1 Overview

Safety Precautions

To avoid risk of electric shock, the following safety precautions must be observed:

- This equipment must only be connected to a supply mains with protective earth
- Do not attempt to service AIRO. Do not attempt to replace the power line cord. There are no user serviceable parts.
- Do not remove equipment covers. Do not attempt to override safety features.
- Do not use AIRO in Oxygen Rich environments or with flammable anesthetics.
- Avoid tripping over or driving over the power and network cables as this can damage them.
- Do not use any electrical cables with damaged insulation or that appear to have been crushed.

AIRO is not waterproof. Water, soap, or other liquids, if allowed to drip into the equipment, can cause electrical short circuits leading to electric shock and fire hazards.

If liquids should accidentally spill into the system, DO NOT apply power or turn the system on until the liquids have dried or evaporated completely.

2.10 Transportation Safety

2.10.1 Overview

Safety Precautions Follow user instructions.

Cautions and Warnings

AIRO's size creates a risk of collision with people and objects. The system is designed to be driven by trained personnel in an average hospital setting with level floors, low thresholds, standard size doors and industrial elevators.

Always drive AIRO with the Pendant Display in the holster. Driving with the display out of the holster creates a danger of both collision and entanglement.

Never allow untrained persons to move or operate AIRO as this creates a hazard of collision and equipment damage. To avoid the hazard of tipping, observe the following safety procedures:

- Never push or attempt to move AIRO from the side. Only push using the hand grips on the system control holster.
- Do not attempt to transport AIRO over damaged flooring.
- Never attempt to transport AIRO with the ring outside of transport position.

Only attempt to transport or operate AIRO on surfaces and inclines designated as appropriate for use by the responsible organization following the criteria established by the manufacturer. Failure to do so creates hazards of collision, tipping, and equipment damage. AIRO is designed to move on level flooring. Inclines cannot exceed 5°. Thresholds must not exceed 9.5 mm (3/8 inch).

Transporting objects, devices or material on the operating tabletop is not permitted as this creates a hazard of equipment damage.

To avoid the potential of serious injury, observe these safety warnings:

- Never transport the system up or down a ramp at an angle. Transporting at an angle on a ramp creates a hazard of tipping which can cause injury.
- Do not attempt to turn while on an incline or decline. Turning creates a hazard of tipping which can cause injury.
- Never unlatch the gimbal while on a ramp. Turning the ring creates a hazard of tipping which can cause injury.
- Never disable the main drive wheel while on a ramp. Disabling the main drive wheel creates a hazard of uncontrolled movement which can cause injury.
- Never store AIRO on a ramp
- Never put AIRO in Scan mode while on a ramp.

Maximum Stopping Distance The maximum stopping distance in a worst-case normal use floor surface is less than 1.0 m.

AIRO weighs over 975 Kg (2150 lbs). Collision creates a risk of injury to personnel or damage to objects.

2.11 Fire Safety

2.11.1 Fire Extinguishers

CO₂ Only

In the event of a fire, only CO₂ fire extinguishers are permitted.

2.12 Power Failure

2.12.1 In Case of Power Failure

Power Interruption AIRO recovers automatically from brief power interruptions (e.g., accidental unplugging). The system shows an error message and reboots when the power is restored.

What Happens if the Power Fails During a Scan

NOTE: You can move the ring away from the table manually.

Power Failure
1. AIRO detects the power outage and aborts the scan.
2. AIRO stores any partial scan.
3. AIRO powers down its battery-driven functions.

How to Restore Limited Functions During a Power Failure



Figure 11

Steps	
1. 	Press Start a on the Control Panel.
2. 	Press On s on the Column keypad or Column remote control.

Functions Available
• The table and column function normally, including putting the table into zero position and releasing it from the column. • AIRO is transportable using internal battery power (same as normal use).

NOTE: Scans cannot be performed until main power is restored.

Battery Failure

All functions cease if the batteries fully discharge. This is highly unlikely.

NOTE: To protect the batteries when unplugged, AIRO powers down automatically when it has been unused for one hour.

Battery Failure
1. Plug AIRO into a wall socket to recharge the batteries.
2. If batteries still fail, contact technical support.

NOTE: Only qualified service personnel may replace the batteries.

2.12.2 Restoring a Session

General Information In the event of a power failure, when power is restored the acquired images will be under Existing Scans.

What Happens When Mains Power is Restored AIRO behaves as if it has been turned on normally.

Options	See
1. Recover data from previous session	6.8
2. Re-scan patient	6 Scanning

3 TRANSPORT MODE

3.1 Chapter Overview

3.1.1 Contents

Topics Covered

Section	See
Controls	3.2
Transportation	3.3
Powering Down	3.4

3.2 Controls

3.2.1 Pendant and Power

Diagram of Controls


Figure 12

No.	Component	Function/Usage
a	Pendant	All controls for transport and scanning except power <i>NOTE: The column has separate controls.</i>
s	Holster	Holds Pendant
d	Hand grips	Molded into holster; user places hands here when driving and to turn
f	Power Controls	Key switch for mains power, Start button to enable pendant
g	USB Port	Data transfer
h	Motion Controls	Move system forward and back (left and right are disabled in transport mode)

3.2.2 Power Switches

Mains Power


Figure 13

No.	Component	Function/Usage	
a	Holster	Holds Pendant	
s	Key Switch		Standby: No functions available except battery charging
			On: Turns on system power
d	Power Indicator Light	Illuminates when the key switch is in the On position and Start g is pressed	
f	USB port	Data transfer	
g	Start	Enables the pendant	
h	Charging Indicator Light	Illuminates when system is plugged to show that the battery is charging.	

3.2.3 Power and Network Cables

Location of Cables The power, network, and ethernet cables attach to the end of the base that is opposite the column.

NOTE: Detach cords before attempting transport.


Figure 14

No.	Components
a	Equipotential stud used to provide a means of linking the earth ground of the system to other equipment to minimize any electrical potential difference between them.
s	Ethernet connectors for View Station (left) and hospital network (right).
d	Power cable

3.3 Transportation

3.3.1 Turn On Power

How to Turn On the System

NOTE: Always fully charge the system before you use it. The best way to ensure this is to leave AIRO plugged in whenever it is stationary.

NOTE: To preserve the batteries, the system automatically shuts off in one hour if left unplugged and unused.



Figure 15

Steps	
1.	Unplug power and network cable(s). The charging indicator light goes off.
2.	 Insert key and turn to on a position.
3.	 Press Start d . • The power indicator light s comes on. • The Pendant displays the Log In screen.
4.	[Option:] Log in on the Pendant.
5.	 Select Transport Mode on the Pendant.

3.3.2 Starting from Parked Position

How to Confirm
AIRO is Ready for
Transport


Figure 16

Steps
1. Power and network cables are unplugged
2. The system is on its wheels
3. Ring is latched in two orientations: • parallel to the base • toward the column

Optional Conditions	Steps	See
System is plugged in	System prompts you to unplug cord(s)	
Ring is in Scan position	1. Push ring away from column. 2. Unlatch. 3. Push to turn the ring until they latch parallel to the base. 4. Push ring toward the column until they latch.	3.3.3
Ring is not latched in one or both orientations	1. Latch ring in Scan position. <i>NOTE: This is to keep the ring from turning while you push it.</i> 2. Push ring away from column. 3. Push to turn the ring until they latch parallel to the base. 4. Push the ring toward the column until they latch.	3.3.4
Base is on the floor, wheels stowed	System prompts you to raise the system.	3.3.5

Universal Option
 Press Home at any time to return to the Choose Mode screen.

3.3.3 Starting from Scan Mode

Using the Latch Use the latch to unlock the ring.


Figure 17

NOTE: The latch locks when the ring is pushed fully into the scan or transport position.

Steps	
1.	Click Transport Mode from the Home screen on the Pendant.
2.	Translate ring along Z-axis away from column to make room to rotate ring for Transport mode.
3.	Lift the latch a to release the lock and rotate the ring to the transport position.
4.	Push ring towards column until it latches.
5.	Unplug power and network cables.
6.	Raise the system to transport. Refer to Section 3.3.5

3.3.4 Rotating and Latching the Ring

How to Rotate and
Latch the ring


Figure 18

Steps	
Rotate the ring until it latches into position:	
1.	• Transport position: Parallel with the base • Scan position: At a 90° angle to the base
2.	(Transport mode only:) Push the ring toward the column until it latches.

3.3.5 Raising the System for Transport

How to Raise the System

Figure 19

Step
Press and hold Up a until the camera window appears. The system is ready for transport. <i>NOTE: System will not transport (camera will not appear) if power cord is plugged in.</i>

3.3.6 Transportation

Diagram


Figure 20

NOTE: Leave the pendant in the holster when transporting the system.

Steps	
1.	Use the display a to see in front of the system.
2.	 Use the buttons s to drive in forward or reverse. NOTE: The right and left buttons are not enabled in transport mode. To turn, keep holding the hand grips on the holster and push right or left.

NOTE: The front-view camera display is only a visual aid to driving. Look around the sides of the system.

Avoid Collisions Transport this equipment following all safety guidelines.

Hand Position Always drive AIRO holding the hand grips on the sides of the holster.

Stopping To stop at any time, release the direction key.

Brakes AIRO engages its brakes when stopped.

AIRO is designed with a fail-safe brake to control movement. In the worst-case on a normal use floor surface, AIRO will stop in less than 1.0 m.

Turning AIRO turns when you manually push left or right on the hand grips. It pivots on its main drive wheel.
NOTE: If you have any doubts about making a turn, stop and evaluate the area first.

Inclines/Declines

The only acceptable position for the system on a ramp is the transport position and parallel with the direction of the incline or decline.

The ring must remain latched in transport position. Inclines and declines must be no more than 5°.

Never unlatch the ring on a ramp.

Never disable the main drive wheel on a ramp.

**Crossing
Thresholds**

AIRO is designed to pass safely over hospital doorway thresholds and into elevators.
It will not cross barriers higher than 9.5 mm (3/8").

Using Elevators

When using elevators:

- Always travel in the elevator with the system.
- Transport only yourself and the system. There can be no additional passengers or equipment.
- Elevator capacity should be at least 1588 kg (3500 lbs).
- Elevator size must be at least 203 cm wide by 165 cm deep (6'8" by 5'5").

3.3.7 After Transport

Options

Options	See
Setup for Scan	3.3.8
Temporary Storage (Park)	3.3.12

3.3.8 Setup for Scan

Determine Location

Where to Place the System

System setup depends on the type of scan being taken or the procedure being performed.

Regardless of the intended procedure, always be aware of the need for audio and visual communication between operator and patient.

- Place the system so the operator and patient can always see and hear each other.
- Always place the system so that it can be easily unplugged in case of an emergency.
- The system and its controls should not restrict the movement of the medical team.
- The ring must have a clear path for movement. Evaluate the area above the scanner as well as the sides.
- The system ventilation system must not be blocked.

Example: Setup During Scan

This diagram is one of many possible setups. The position of **Airo**, accessories, medical personnel, and the patient varies according to need.


Figure 21

Components

No.	Component
a	CT technician holding pendant
s	Protective barriers
d	Patient on operating table
f	Gimbal
g	Ring
h	Other medical personnel

How to Ensure Ventilation

Air must be allowed to circulate through the vents to ensure proper operation and to avoid overheating.

Do not block or cover ventilation openings on the system, for example, with drapes. This creates a fire hazard.


Figure 22




Step
Allow air to circulate freely around the intake and exhaust vents a in scan mode.

3.3.9 Lowering the System

How to Lower the System

Figure 23

Steps	
1.	 Select Menu on the Pendant Transport Screen to bring up the Main Menu options.
2.	 Select Scan on the Main Menu .
3.	 Press Down until the base is on the floor.
4.	Plug in the power cord.
5.	If applicable, plug in the network cable.

Be sure there are no objects under the base when it is being lowered, such as power cords.

3.3.10 Unlock and Rotate the Ring

How to Unlock and Rotate the Ring


Figure 24

Step
Pull the release latch to unlock the ring and rotate it manually. The ring will lock in at $\pm 90^\circ$. <i>NOTE: Rotate the ring so that the window inside the ring faces the forward direction of the scan. The scan will either be conducted toward the column or away from the column.</i>

3.3.11 Setting the Orientation for Imaging

How to Set the Orientation for Imaging

The ring rotates $\pm 90^\circ$. This allows users to scan toward either end of the patient support.

Step

Turn the ring manually so that the window inside the ring faces the forward direction of the scan.

Next Step

Step	See
Transfer patient using Column, Table, and Shuttle	4 COLUMN AND TABLE 5 TRANSFERRING THE TABLE

3.3.12 Temporary Storage

How to Park

Leave AIRO on its wheels when you are not scanning. Only lower the base to the floor for scanning.

NOTE: Park on level surfaces only. Never store on a ramp.


Figure 25

Steps	
1.	 Select Home to exit camera view.
2.	 Turn the key switch to the standby position.
3.	Plug in the power cord.

3.4 Powering Down

3.4.1 Set to Standby Mode

Overview

Whenever AIRO is not in use it should be in Transport position, plugged in, and set to **Standby**.

How to Enter Standby Mode

Steps	
1.	Leave the system on its wheels.
2.	Turn the ring so it is aligned with the base.
3.	Push the ring toward the column until it locks into place.
4.	Plug in the power cord. The charging indicator light illuminates.
5.	Insert the key and turn it to Standby. The indicator light blinks while the PC shuts down.

NOTE: Always keep AIRO plugged in to maintain the batteries.

4 COLUMN AND TABLE

System AIRO could be connect with surgical tables
from TRUMPF Medizin Systeme GmbH +Co. KG

4.1 Chapter Overview

4.1.1 Contents

Topics Covered

Section	See
Precautions	4.2
Overview of the Column	4.3
Overview of the Controls	4.4
Column Keypad	4.5
Wireless Remote Control	4.6
Controlling the Tabletop	4.7
Functions on Wireless Remote Control Only	4.8

4.2 Precautions

4.2.1 Overview

Before You Begin

In addition to the safety procedures described in 2 *SAFETY*:

Secure the patient before adjusting the operating tabletop and provide active support during adjustments.

Only execute operating tabletop functions when there is no shuttle against the operating tabletop column.

Risk of patient injury from unintended movement. There is no Emergency Stop function on the Column Control or on the wireless Remote Control.

NOTE: The Emergency Stop is on the pendant. See 2.3 In Case of Emergency.

4.3 Overview of the Column

4.3.1 Surgical Column

Description	The surgical column for system AIRO supports the operating tabletop. The column is used with these accessories: • Wireless Remote Control and Charger • Operating tabletop (various models) • Tabletop segments and accessories • Patient shuttle (various models) <i>NOTE: For more information on accessories, see 11 ALLOWEDEQUIPMENT.</i>

4.3.2 Column

Surgical Column


Figure 26

No.	Component
a	Point of electrical contact to operating tabletop
s	Holding wedge for operating tabletop
d	Column keypad
f	Lowest column cover (column base)
g	Base plate

4.3.3 Audio Signals

How to Interpret Audio Signals

The **Airo** surgical column emits audio signals. A corresponding message appears on the wireless Remote Control display.

Action or Status	Sound
Turning on	A series of ascending tones
Turning off	A series of descending tones
Arriving at selected position	A single tone
Tone confirming an action	Two tones
Warning or error	Three high-pitched tones
Emergency stop	Series of descending tones

4.4 Overview of the Controls

4.4.1 Column Keypad and Wireless Remote Control

What Controls Can Be Used These control units can be used to adjust the operating tabletop:

Control Unit	See
Column Keypad	4.5
Wireless Remote Control	4.6

NOTE: For step-by-step instructions on performing specific functions, see 4.7 Controlling the Tabletop.

4.4.2 Using the Controls

How to Use the Controls *NOTE: Some functions stop automatically. Functions such as **Zero Position** stop when the position is reached.*

Each key function is identified by a corresponding pictogram.

Option	Action
To adjust the tabletop or column	Press the function key(s) until the desired position is reached
To stop the movement	Release the key
To continue a function beyond an automatic stop	Briefly release the function key and press it again

Using the Column Without the Tabletop When there is no tabletop on the column, only **Zero Position**, **Lift**, and **Lower** function are operable.

4.5 Column Keypad

4.5.1 Overview

Diagram


Figure 27

Most Commonly Used Functions

No.	Function	See
a	On	4.7.3
s	Lock and Unlock Two-Key Function <i>NOTE: The column keypad requires you to press two buttons to perform a function unless you specifically unlock this safety feature. See 4.7.4 Lock and Unlock Two-Button Operation on the Column.</i>	4.7.4
d	Raise/lower column	4.7.7
f	Brake for column rotation	4.7.5
g	Zero Position	4.7.8

Tabletop Positioning Functions


Figure 28

No.	Tabletop Movements	See
a	Trendelenberg/Anti-Trendelenberg	4.7.11
s	Tilt left/right	4.7.10
d	Raise/lower leg section	4.7.12
f	Longitudinal displacement at the head/feet	4.7.12
g	Raise/lower lower back section	4.7.12
h	Raise/lower upper back section	4.7.12

Other Functions


Figure 29

No.	Function
a	Set head position
s	Malfunction
d	Activate emergency mode: See 12.6.3 <i>Table Will Not Release</i>
f	Locate wireless Remote Control

4.5.2 Column Keypad Visual Signals

Key to Visual Signals

State		Meaning
	Illuminated	On
	Illuminated	Head position is left
	Illuminated	Head position is right
	Illuminated	Tabletop is at zero position
	Blinking	Emergency state is on
	Steady green: Blinking green:	Brake is on. Tabletop cannot be rotated. Brake is off. Tabletop can be rotated.
	Illuminated	Malfunction Use is limited Contact Customer Service

4.6 Wireless Remote Control

4.6.1 Overview

**Remote Control
Detects the System**

The wireless Remote Control automatically detects the AIRO system:

- When the column is on (see *4.7.3 On and Off*) and
- The Remote Control is within range of its signal

The wireless Remote Control works with only one AIRO system at a time and must be paired with that system. See *4.8.7 Service* or contact technical support.

4.6.2 Charging the Wireless Remote Control

How to Charge the
Wireless Remote
Control


Figure 30

Step
Place the wireless Remote Control in the charging cradle as shown.

4.6.3 Overview of the Wireless Remote Control

Diagram The wireless Remote Control includes a status bar, touchscreen, and keypad.



Figure 31

No.	Area	See
a	Touchscreen	4.6.4
s	Column controls	4.5
d	Tabletop movement controls	4.7
f	Zero Position, brake, and off	4.7

4.6.4 Touchscreen

Status Bar Key The top line of the touchscreen is a status bar.

The status bar displays the following information:

Icon	Refers to	State	Meaning
	Battery Charge	White bar:	The length of the white bar shows the level of charge.
		White moving bar:	Battery is charging
		Blinking red:	Battery is empty. Place Remote Control in charger.
	Tabletop Position	White:	"Zero" (starting) position
		Orange:	Starting position error. Contact Customer Service.
	Brake	White:	Brake is on
		Blinking orange:	Brake is off
	Wireless Signal	White:	Signal contact OK
			Signal is blocked
		Orange:	Remove obstacles between tabletop and Remote Control.
		Absent:	No signal. • Tabletop is out of range • Tabletop is turned off • Connection is not established between Remote Control and tabletop. See 4.8.7 Service.

Touchscreen Menu Key


Figure 32

Option	Function	See
Slow	Adjust the operating tabletop at reduced speed.	4.8.2
Position	Select tabletop position or define a new one.	4.8.5
Service	Change settings (e.g., language).	4.8.7
Flex	Simultaneously execute Anti-Trendelenburg and Lower lower back section functions.	4.8.3
Reflex	Simultaneously execute Trendelenburg and Raise lower back section functions.	4.8.3
Light	Turn on the light under the tabletop.	4.8.4

Touchscreen Selectors

Use the two buttons below the touchscreen to cancel or confirm touchscreen selections.

Select	To
	Cancel or go back
	Confirm selection

4.6.5 Keypad

Diagram of the Keypad

The keypad functions change with the type of operating tabletop attached. Only the functions that are achievable with each tabletop model are illuminated.


Figure 33.

No.	Key	See
a	Raise	4.7.7
s	Lower	4.7.7
d	Tabletop adjustments	4.7.5
f	Zero Position	4.7.8
g	Brake	4.7.5
h	Off (turns wireless remote off)	4.7.3

Key to Tabletop Adjustment Functions

Instructions to use these functions are in *4.7 Controlling the Tabletop*.


Figure 34

No.	Tabletop Adjustments
a	Tilt left/right
s	Trendelenburg /anti-Trendelenburg
d	Select joint, left/right
f	Raise/lower leg sections
g	Raise/lower leg section or the motorized seat section extension
h	Transverse displacement, left/right
j	Longitudinal displacement at the head/feet
k	Raise/lower the lower back section
l	Raise/lower the upper back section
;	Raise/lower the motorized head rest

4.7 Controlling the Tabletop

4.7.1 Overview

Common Functions Most functions are available on both the column keypad and the wireless Remote Control.

*NOTE: The exception is the **On** function, which is only on the column keypad.*

See 4.8 for the functions specific to the wireless Remote Control.

4.7.2 Instructions

Where to Find Instructions for Functions

Function	See
On and Off	4.7.3
Lock and Unlock Two-key Operation	4.7.4
Rotating and Braking	4.7.5
Lift and Lower	4.7.7
Zero Position	4.7.8
Head Position	4.7.9
Tilt	4.7.10
Trendelenburg and Anti-Trendelenburg	4.7.11
All other tabletop functions	4.7.12
Remote only:	4.8.1
• Slow	4.8.2
• Flex and Reflex	4.8.3
• Light	4.8.4
• Position	4.8.5
• Stored Positions	4.8.6
• Service	4.8.7

Instructions for Transferring Patients and Tabletops

Instructions to transfer tabletops between the column and the shuttle are in 5 *TRANSFERRING THE TABLE*.

4.7.3 On and Off

How to Turn On the Column Be sure the system is plugged into wall power.:

Column Keypad	Remote Control
 Press On for more than 2 seconds.	There is no equivalent command on the wireless Remote Control. <i>NOTE: The remote only communicates with the column when the column is on.</i>

How to Turn Off the Column Immediately To turn off the column immediately:

Column Keypad	Remote Control
 Simultaneously press On and Off for more than 2 seconds	There is no equivalent command on the wireless Remote Control. <i>NOTE: The Off icon on the wireless remote turns off the remote.</i>

Automatic Shutdown

When the column detects it is not being used, it automatically shuts off.

Automatic Shutoff	Time Since Last Command
Tabletop on column	Four hours
No tabletop on column	One hour

4.7.4 Lock and Unlock Two-Button Operation on the Column

Safety Feature The column keypad requires two keys to perform a function unless the user specifically unlocks it.

Option	Action
Locked <i>NOTE: Default: Two-button operation.</i>	User must press On + function to perform any function.
Unlocked	User can perform functions with one button.

How to Unlock To enable one-key function on the column keypad:


Step	Action
 Press for two seconds.	Unlock illuminates. Automatically locks 15 seconds after last button is pressed.

Instructions
Assume Keypad is Locked *NOTE: The instructions in this manual assume that the column keypad is locked. If you have unlocked it, skip the “press **On**” instruction.*

4.7.5 Rotate and Brake the Column

4.7.6 Instructions

How to Rotate and Brake the Column

Column Keypad	Remote Control
	

Tabletop is Locked	Tabletop is Unlocked
1. Press On and Brake .	Press Brake .
2. Push the tabletop to rotate the operating tabletop column. <i>NOTE: The operating tabletop may be rotated 220° right or left.</i>	
3. Press On and Brake again to lock the column into position.	Press Brake again to lock the column into position.


NOTE: If you do not apply the brake, the brake will automatically lock after 30 seconds.

4.7.7 Lift and Lower

Definition and Cautions

The lift and lower functions raise and lower the entire operating tabletop. The range of adjustment is 550 mm (21.7").

How to Use the Lift Function


How to Use the Lower Function

To avoid equipment damage, move any lowered leg sections to it's respective zero position or level position before lowering the tabletop.

To lower the tabletop, press:

4.7.8 Zero Position (Transfer Position)

Definition and Cautions

Zero position is the starting position for table and column. It is the setting used when transferring tabletops to and from the column.


Pressing **Zero Position** automatically moves all the electrically adjustable operating tabletop positions to their defined starting (neutral) positions.

Monitor the patient's position. When approaching the zero position, inclined patient positions change. This may create a hazard to the patient and danger of collision.

How to Use Zero Position

There is a **Zero Position** key on the column keypad and on the wireless Remote Control. Use only one controller for the function.

Steps	
1.	Adjust any swivel leg sections so that they are not on top of one another.
2.	Position accessories so that collisions can be ruled out when approaching the zero position.
3.	 Press Zero Position. <i>NOTE: The tabletop and column stop automatically when they reach zero position.</i>
Column Keypad	Remote Control
	
4.	Manually place the section segments of the mechanically adjustable joints in the zero position (the horizontal position) (align the gear markings).
5.	<i>Tables with mechanically adjustable longitudinal displacement:</i> Manually move the operating tabletop in the zero position in the longitudinal direction.

4.7.9 Select the Head Position: Column Keypad Only

Airo Detects the Head Position

Airo detects the head position based on the direction of the operating tabletop (head/foot end) on the column. The patient's head position is illuminated on the column keypad.

NOTE: There is no equivalent function on the Remote Control.

Figure 35

NOTE: The head position of the patient must match the display on the column keypad.

If you have placed the patient on the tabletop in the reverse orientation, override the detected head position.

How to Override the Detected Head Position

When the column keypad displays **a** or **s** do not match the patient's head position, adjust the setting for the head position to match the actual patient head position.

Using the keypad on the column:

Options	Instructions	
To set the tabletop head position LEFT		Press LEFT
To set the tabletop head position RIGHT		Press RIGHT

4.7.10 Tilt


Definition and Limits

When tilted, the operating tabletop is inclined to the left or right along its longitudinal axis.


NOTE: For safety reasons, the tilt function automatically stops at 15°.

How to Use Tilt Functions

To Tilt left, press:

To Tilt right, press:

Options	Instructions
Stop at 15°	Automatic
Adjustments beyond 15°	Briefly release the function key and press it again.
Tilt limitations in extreme Anti-Trendelenburg positions	Increase the degree of tilting by reducing the Anti-Trendelenburg position.
Maximum tilt: 30°	

4.7.11 Trendelenburg and Anti-Trendelenburg

Definition and Cautions

With the Trendelenburg functions, the operating tabletop is tilted in the head or foot direction along its transverse axis.

Range of Adjustment
$\pm 45^\circ$
$\pm 65^\circ$ for operating tabletops with an additional Trendelenburg adjustment

May create a hazard to the patient and danger of collision. Clear the area beneath the tabletop before executing function.

If necessary, move the leg sections to zero position before executing the function.


How to Use the Trendelenburg Function

To move into the Trendelenburg position, press:

Column Keypad	Remote Control

How to Use the Anti-Trendelenburg Function

To move into the Anti-Trendelenburg position, press:

Column Keypad	Remote Control
	

4.7.12 Other Tabletop Functions

Other Functions Available

All other table functions depend on the tabletop and accessories that are installed.

Read the user manuals specific to the tabletops and accessories for descriptions of their functions.

The available functions are illuminated.

How to Use Other Functions

Option	Action
Adjust the tabletop or column	Column Keypad: Press On + function key until the desired position is reached Remote Control: Press the function key until the desired position is reached
Stop	Release the key
Continue a function beyond an automatic stop	Briefly release the function key and press it again

4.8 Functions on Wireless Remote Control Only

4.8.1 Functions Specific to Wireless Remote Control

Functions Only on
Wireless Remote
Control

Option	See
Slow	4.8.2
Flex and Reflex	4.8.3
Light	4.8.4
Position	4.8.5
Stored Positions	4.8.6
Service	4.8.7

4.8.2 Slow

How to Use the Slow Function

The **Slow** function reduces all mechanized tabletop movements by 50%.


Figure 36


Step
Touch Slow to turn the function on or off.

4.8.3 Flex and Reflex

Definition The **Flex** and **Reflex** functions move two tabletop sections simultaneously to create a bend between the seat and back sections of the operating tabletop.



Options	Movements
Flex	Lowers the lower back plate, Anti-Trendelenburg
Reflex	Raises the back plate, Trendelenburg

How to Use Flex and Reflex Use the touchscreen menu:


4.8.4 Light

How to Use the Light Function

Steps	
1.	 Press Light.
2.	 Endolight appears.
Options	
To light the area below the tabletop	1.Press Endolight 2.Press OK 3.You return to the Main Menu.
To turn the light off	From the Main Menu: 1.Press Light 2.Press Endolight . 3.Press OK You return to the Main Menu.
To exit	From Endolight : 1.Press Back to return to the Main Menu

4.8.5 Position

How to Use the Position Menus

Steps	
1.	 Press Position.
2.	 Scroll up or down the submenu with the arrows. The menu item is highlighted.
Options	
Press OK to select the highlighted item.	
Press Back to exit the submenu and return to the Main menu.	

**Position Submenu
Options**

Options	Step
To save the current position of the table	Select Save Current Position. See this section for directions.
To select a stored position	Select a category and refer to 4.8.6 Stored Positions : • Main • General/Visceral: Internal Organs • Ortho • Trauma • Gyn/Uro • Neuro

NOTE: A number in parentheses is the number of stored positions in that category.

**How to Save the
Current Position**

Current positions are always saved in the **Main** category. See **4.8.6 Stored Positions** to change the category of a stored position

NOTE: Positions using two additional motorized components with operating tabletop U26/ST26 cannot be stored.

Steps
1. Adjust the operating table to the desired position.
2. Touch Position .
3. Scroll to Save Current Position .
4. Press OK .

The position is assigned a consecutive number (e.g., P_001) and is stored in the **Main** category.

To move the position into another category, see **4.8.6 Stored Positions**.

**Positions are
Automatically
Numbered**

The system assigns the position a number.

Stored positions can be renamed by Technical Service or the responsible on-site hospital technician.

**How to Select a
Stored Position**

Steps
1. Touch Position .
2. Use the arrows to scroll to a category. Press OK .
3. Use the arrows to scroll to a position number. Press OK .
4. Select Move to Position .
5. Press and hold OK until a signal tone sounds and Target Position Reached appears on the display.

NOTE: You can only select positions that use the installed table and accessories.

4.8.6 Stored Positions

General Up to 30 positions can be saved in predefined categories.

How to Change the Category of a Position



Steps	
1.	Touch Position .
2.	Use the arrows to scroll to a category. Press OK .
3.	Use the arrows to scroll to a position name. Press OK .
4.	Scroll to Change Category . Press OK .
5.	The current category is marked with an asterisk(*). Use the arrows to scroll to the new category.
6.	Press OK . Category Changed appears on the display.

How to Delete a Stored Position

Steps	
1.	Touch Position .
2.	Use the arrows to scroll to a category. Press OK .
3.	Use the arrows to scroll to a position name. Press OK .
4.	Select Delete Position .
5.	Press OK to permanently delete the stored position. Press Back to keep the position in memory and go back.

4.8.7 Service

How to Use the Service Menu

Steps	
1.	 Press Service.
2.	 Scroll up or down the submenu with the arrows. The menu item is highlighted.

Options

Press **OK** to choose the highlighted item.

Press **Back** to exit the submenu and return to the higher menu.

Service Submenu Options

Task	Select
To change the language	Language
To link the remote with a table	Operating Table
To set or release the shipping lock	Shipping Lock
To view messages	Recent Messages
To remove all stored positions	Delete Entire Position Memory
To temporarily slow table movements	Temporary Slow Traverse
To access system functions	System (password-protected)

How to Set the Language

Steps	
1.	Touch Language .
2.	The current language is marked with an asterisk (*). Use the arrows to scroll to the new language.
3.	Press OK to change the language. Press Back to exit the submenu and return to the higher menu.

How to Link the Remote with an AIRO System

For safety reasons, the remote control can be used with only one system at a time.

To link the remote with a system:

Steps
1. Follow the steps in 6 SCANNING to turn on the system and enter Scan Mode.
2. Follow the steps in 4.7.3 On and Off to turn on the column.
3. On the remote, touch Operating Table .
4. The display shows the serial number of the operating column to which it is linked. (The serial number of the column is on the rating plate. See 9.16.1 Overview) If None is displayed, the Remote Control needs to be linked.
5. Press OK to list the systems that can be accessed by the remote control.
6. Use the arrows to scroll to the serial number.
7. Press OK to confirm. Press Back to return to the higher menu.

The display shows the wireless signal. The keys for available functions illuminate.

How to Unlock the Remote

The remote control is delivered locked.

Step
Touch Shipping Lock .

If the remote control is locked, the display says, **Turn off Shipping Lock?**

If the remote control is not locked, the display says, **Activate Lock?**

Options
To release the shipping lock, press OFF . Press OK to confirm.
To lock the remote, press OK to set the lock.
To leave the setting as it is, press Back to return to the Service menu.

How to View Recent Messages

Steps
1. Touch Recent Messages .
2. The display shows the last ten error and status messages. The most recent message is first.
3. Use the arrows to scroll through the messages.
4. Press Back to return to the higher menu.

How to Clear the Entire Position Memory

Steps
1. Touch Delete Entire Position Memory .
2. Press OK to delete every stored position from memory. Press Back to return to the higher menu.

How to Set a Temporary Slow Traverse

Use this menu item to temporarily slow the speed of the operating table functions to 50% or 75% of the normal speed for up to 90 seconds.

Steps	
1.	Touch Temporary Slow Traverse. The current slow speed is marked with an asterisk (*)
2.	Use the arrows to scroll to a speed. Press OK to select the slow speed setting. Press Back to exit the submenu without changing the speed.
3.	Use the arrows to scroll to a deactivation time. The current deactivation time is marked with an asterisk (*) Press OK to select the deactivation time. Press Back to exit the submenu and return to the higher menu.

System Submenu

This menu item is password-protected and accessible only by authorized service technicians.

5 TRANSFERRING THE TABLE

5.1 Chapter Overview

5.1.1 Contents

Topics Covered

Section	See
Transfer Safety	5.2
Transferring to or from the Shuttle	5.3

5.2 Transfer Safety

5.2.1 Introduction

Transferring Patients

AIRO allows different patient supports to be used with the system, but most often the patient will be on a surgical table.

This chapter describes transferring the surgical table from a shuttle to the column. It also describes how to remove the table from the column and put it back on the shuttle.

Before using the table or other patient support, or the shuttle, read the applicable instruction manual.

5.2.2 Patient Transfer Safety

General Precautions	Monitor the transfer carefully. The stipulated sequence of actions has to be followed, checked and carried out with care.

Limited Functions	For safety reasons, the operating table movements in the Transfer function are limited. You can only use them: • When the table is in the zero position • When the sensor monitoring on the column authorizes the transfer

Risk of table tipping. Risk of injury to the patient and damage to equipment.

NOTE: Remove or fold in the accessories before transporting the operating table or tabletop. Align the spreadable leg sections parallel to the longitudinal direction of the operating table.

Precaution on Patient Position	Take care to position the patient in a central position on the table. The center of gravity must be over the seating wedges.

NOTE: The transfer cannot take place if the weight is not centered.

Risk of injury to the patient. Extreme shifting of the center of gravity during the transfer may prevent the column's holding wedges from locking into the tabletop.

5.3 Transferring to or from the Shuttle

5.3.1 Transferring the Table from the Shuttle to the Column

Only Use One Control The transfer procedure is carried out using the **Transfer** pedal on the shuttle or using the **Zero position** key on the remote control.

NOTE: Do not use the pedal and remote control at the same time.

How to Stop the Transfer To stop the process, release the pedal or the key.

How to Restart a Stopped Transfer Restart the process by re-actuating the **Transfer** pedal on the shuttle or the **Zero position** key on the remote control.

How to Transfer the Table to the Column To transfer the table to the surgical column:

Figure 37


Figure 38

Steps	
	Set the shuttle to the zero position.
1.	<i>NOTE: When using shuttle types with angle and height adjustment, transfer is only possible when the holding wedges on the shuttle are horizontal and the shuttle column is at its highest position.</i>
2.	Set 5th wheel pedal a to Freewheel.
3.	Set the Brake pedal S to Unbraked.
4.	Put the shuttle against the column.
	 Risk of material damage due to collisions. When running the shuttle up to the operating table column, the operating tabletop may collide with the column. If necessary, bring any lowered section segments on the operating tabletop to a horizontal position and raise or lower the column accordingly. For more information, consult the instruction manual for the operating tabletop you are currently using. <i>NOTE: When docking on the column, always ensure that the wheel is locked in behind the column.</i> Shuttle present, Transfer mode appears on the remote control display.
5.	Press and hold down the d Transfer pedal on the shuttle or the Zero Position key on the remote control.
	The column moves to the upper lift position.
6.	Transfer the operating table to the column.
	<i>NOTE: Ensure that the operating table or other patient support is properly locked onto the holding wedges of the column.</i> A signal (two tones) indicates the end of the transfer and Tabletop (type) recognized appears on the remote control display.
7.	Move the shuttle away after the signal sounds.

Error: Table is Not Locked

Step
Raise the operating tabletop on the “hanging” side until the wedge holders audibly lock into place

5.3.2 Transferring the Table to the Shuttle

How to Transfer the Table to the Shuttle To transfer the table from the column back to the shuttle:

Figure 39


Figure 40

Steps	
	Set the shuttle to the zero position
1.	<i>NOTE: When using shuttle types with angle and height adjustment, transfer is only possible when the holding wedges on the shuttle are horizontal and the shuttle column is at its highest position.</i>
2.	Set a 5th wheel pedal to Freewheel.
3.	Set the s brake pedal to Unbraked.


Figure 41

4.	Put the shuttle against the column. <i>NOTE: When docking on the column, always ensure that the wheel is locked in behind the column.</i> The remote control displays Shuttle present, Transfer mode.
5.	Press and hold down d the Transfer pedal on the shuttle or the Zero Position key on the remote control. The column moves to the zero position and then to the lower transfer position, lowering the tabletop onto the shuttle.
6.	Ensure that the table is properly locked onto the holding wedges of the column. A signal (two tones) indicates the end of the transfer. The remote control displays Tabletop removed.
7.	Move the shuttle away after the signal sounds.

5.3.3 Reversing Direction

How to Reverse the Direction During a Transfer To interrupt a transfer and reverse direction:


Figure 42

Steps	
1.	To interrupt the transfer, release the Transfer pedal d or Transfer key on the keypad. The remote control shows the message Transfer interrupted, change direction?
2.	To reverse movement, confirm by pressing the OK key or by briefly pressing any key on the column keypad. <i>NOTE: Do not use the level and the keypad at the same time.</i> The change in direction of travel is confirmed by a signal tone (single tone). The remote control displays the message Direction of transfer changed.
3.	Press the Transfer pedal d on the shuttle or the Zero Position key on the remote control again. The movement takes place in the opposite direction.

6 SCANNING

6.1 Chapter Overview

6.1.1 Contents

Topics Covered

Section	See
System and Controls	6.2
Before Scanning	6.3
Moving the Ring in Scan Mode	6.4
Prepare Scan Field of View	6.5
Welcome Menu	6.6
Typical Scan Sequence: New Scan	6.7
View an Existing Scan	6.8
Emergency Scanning	6.9
After Scanning	6.11
Service Menus	6.12

6.2 System and Controls

6.2.1 AIRO System

Diagram of System
Configured for
Scanning

Figure 43

No.	Component	Function/Usage
a	Ring	Scans patient and sends data to pendant
s	Window	Opening through which scan is taken
d	Base	Supports and moves the system, includes power and data connections. <i>NOTE: Casters are retracted into the base.</i>
f	Gimbal	Supports the ring
g	Latch	Releases gimbal, enabling rotation
h	Column	Supports some models of surgical tables/chairs (see parts list for specific part numbers)
j	Pendant Display	Controls system movement, imaging, and data storage <i>NOTE: The Pendant Display is shown stored in the holster.</i>

Details

In scanning mode, the base is lowered and locked. The ring is around the surgical table. The system operator has the Pendant Display out of its holster and stands back from the system and/or behind a protective barrier.

The operator runs the system and views the images on the Pendant Display.

During a scan, the ring moves down the patient axis, then returns to the end position to upload the image data.

6.2.2 Pendant Display

Diagram ofPendant
Display


Figure 44

No.	Component	Function/Usage
a	Touchscreen	Touch to enter data or select functions Serves as display for camera view
s	Emergency Stop	Press to halt all scanning and mechanical motion functions
d	Scan	Press to conduct a scan
f	Directional Keys	Used to move ring and gimbal

Details

The CT operator uses the Pendant Display to move the ring and gimbal, turn the laser pointer on and off, set the parameters for each scan, conduct the scan, and interact with the resulting images.

Options	See
Move the ring and gimbal	6.4
Scan	6.7
Use Laser for Positioning	6.7.11
Move the table and column	4

6.3 Before Scanning

6.3.1 Overview

Before Scanning

Steps	See
1. Follow the instructions to set up for scan	3.3.8
2. Confirm ring is facing the direction of the intended scan before moving the patient to the column	6.7
3. Transfer patient and patient support to the surgical column	5
4. Option: Log in	6.3.2
5. Choose Scan mode	6.3.3
6. Prepare scan field of view	6.5

6.3.2 Option: Logging In

How to Log In

Your system may be set to time out or log out for security. In this case:


Figure 45

Steps	
1.	Enter password.
2.	Press Next .

6.3.3 Choosing Scan Mode

How to Choose Scan Mode


Figure 46

Step
Press Scan a on the touchscreen.
<i>NOTE: If the system is not lowered, you are prompted to lower the system. See 3.3.9.</i>

Required Maintenance Scans If the system detects that Required Maintenance Scans have not been performed recently, you are prompted with a dialog box indicating which Scans need to be performed.


Figure 47

Condition	Options
Warm-up Scan	A Warm-up scan must be performed before any other scans can be performed. <i>See 6.12 Service Functions</i>
E-stop test	E-stop test must be performed to ensure scanner safety overrides are functioning properly in case of emergency. <i>See 6.12 Service Functions</i>
Gain Calibration	You must calibrate the system to ensure the best image quality possible. <i>See 7.2 Daily Gain Calibration.</i>
Quality Check Scan (SAB)	Image Quality testing must be performed per hospital requirements. (SAB) indicates Scan Away from Base. <i>See 6.12 Service Functions.</i>
Quality Check Scan (STB)	Image Quality testing must be performed per hospital requirements. (STB) indicates Scan Towards Base. <i>See 6.12 Service Functions</i>

6.3.4 Connecting System

How to Connect to Power and Hospital Network


Figure 48

Steps	
a	Connect power cord to standard wall outlet.
s	Option: Connect ethernet wire to hospital network.
d	Option: Connect to surgical navigation system.
f	If required for electrostatic protection from other equipment, connect equipotential stud f to ground.
Press Next .	

Task Bar icons	
g	Ring orientation
h	Power
j	Hospital ethernet connection

Task Bar icons	
K	Navigation ethernet connection
	System temperature
	<i>NOTE: System temperature icon indicates status of x-ray tube.</i>
	Cyan - tube is cold, a warm-up scan must be performed
	Green - 0 to 33% of tube heat capacity is used
	Yellow - 33% to 66% of tube heat capacity is used
	Orange - 66% to 100% of tube heat capacity is used
	Red - 100% of tube heat capacity is used, a cool down wait time is enforced

6.4 Moving the Ring in Scan Mode

6.4.1 Translating the Ring (z axis)

How to Translate the Ring on the Base

When in Scan mode, use the direction keys on the Pendant Display control ring translation.


Figure 49

Step	Action
To move the ring back and forth	Press and hold the left a or right s directional button on the Pendant Display

6.5 Prepare Scan Field of View

6.5.1 Position and Diameter of Scan Field of View

Scan Field of View


Figure 50

The scan field of view is a 51.2 cm (20.2") diameter circle at the center of and parallel to the bore.

6.5.2 Remove Non-Essential Objects

Non-Essential Objects

Remove all non-essential objects and equipment from the ring bore and the scan field of view before scanning. When necessary, place IV tubes where the motion of the system will not affect them.

Extra objects in the ring bore, even when they are out of the scan field of view, may degrade the scan image.

6.5.3 Remove Metal Objects

Keep Metal Objects Away from System Be sure there are no metal or other radiopaque objects in the bore or near the scan field of view.
NOTE: This includes table accessories such as arm boards.

Metal objects in or near the scan field of view can degrade the scan image.

6.6 Welcome Menu


6.6.1 Welcome Menu and Functions

Options on the Welcome Menu

Figure 51

Step		
Select a task by touching the icon.		

Option	Function	See
New Scan 	Select to start the scan workflow.	6.7.1
Existing Scans 	DICOM standard viewer to display existing images that are on the system	6.7.3
Service Functions 	Service Functions Menu	6.12

Option	Function	See
 Quick Scan	Perform Emergency Scan (skip patient data entry)	6.9
 Menu	Opens additional options menu	6.6.2

6.6.2 Menu

Universal Option:
Menu


Figure 52

Step
Press Menu a at any time to access the Menu.
Options
Select s Home to return to the Choose Mode screen. See 6.3.3.
Select d New Patient to enter information for a new patient. See 6.7.
Select f Transport to exit scan mode and enable system transport.
Select g Warm Up to initiate warm up scan sequence.
Select h Info to learn more about AIRO.
Select j System Lock to enable security requiring a password to access system functions. See 6.3.2

6.7 Typical Scan Sequence: New Scan

6.7.1 New Scan

How to Initiate New Scans



Figure 53

Step	
	Select New Scan a. AIRO displays the Patient Data screen.

Select Source of Patient Data







Figure 54

Options	See
	Select New Patient a to create a new patient record on AIRO. 6.7.2
	Select Existing Patients s to access the list of patient records currently on AIRO. 6.7.3
	Select Hospital d to access the hospital network: • Select a patient record and download basic information • View the work list 6.7.4
	Select Navigation f to retrieve patient information from the navigation system. 6.7.5

6.7.2 Entering New Patient Data

How to Enter New Patient Data

Figure 55


Steps	
1.	 Press New Patient from the Patient Data Screen.
2.	Enter name and date of birth; select gender a .
3.	Press Next . See 6.7.6 <i>Confirming Ring Rotation</i> .


How to Correct Errors

Problem	Correction
Error in patient name	Edit the name using the edit function under Existing Patients. See 6.7.3

6.7.3 Select Existing Patient Data

How to Select a Patient With Existing Scans

Steps	
1.	 Press Existing Patient from the Patient Data Screen See 6.7.1.
2.	Select an existing patient name from the list a or search by patient name or ID S .
3.	Editor delete patient information if required d .
4.	Press Next . See 6.7.6 <i>Confirming Ring Rotation</i> .

6.7.4 Entering Patient Data from Hospital Worklist

How to Enter
Patient Data from
Hospital Worklist


Figure 56



Figure 57

Steps	
1.	 Press Hospital from Patient Data Screen See 6.7.1.
2.	Select a a server and a worklist.
3.	Select a patient name from the list or search by patient name or ID.
4.	Enter the patient's weight.
5.	Press Next . See 6.7.6 <i>Confirming Ring Rotation</i> .

6.7.5 Entering Patient Data from Navigation System

How to Enter Patient Data from Navigation System


Figure 58

Steps	
1.	 Press Navigation from the Patient Data Screen See 6.7.1.
2.	Select a a navigation system and a worklist.
3.	Select a patient name from the list or search by patient name or ID.
4.	Enter the patient's weight.
5.	Press Next . See 6.7.6 <i>Confirming Ring Rotation</i> .

6.7.6 Confirming Ring Rotation

How to Confirm the Ring Rotation

Figure 59

Step
Confirm the ring rotation for imaging. <i>NOTE: The window in the bore faces the direction of the scan.</i>
Options
If the ring is not locked in the $+90^{\circ}/-90^{\circ}$ position, you are prompted to lock the ring.
If the ring is in the wrong orientation, move the table and gantry and reposition the ring.
To change parameters, press Back .
To proceed, press Next .

6.7.7 Selecting Exam Region

How to Select the Exam Region


Figure 60

Steps	
1.	Select main body region to be scanned. (left column)
2.	Select sub-region to be scanned. (right column)
3.	Adjust Kernel and Field of View if required.
4.	Select with or without Navigation. (if available)
5.	The display shows the kV and mA for the proposed scan.
Options	
To change the ring orientation, press Back .	
To proceed, press Next .	

6.7.8 Selecting Navigated Scan

How to Select a Navigated Scan

Figure 61

Step
Select Navigation.
<i>NOTE: Navigation system must be connected and configured properly when utilizing the Navigated scan option.</i>
Options
To change parameters, press Back .
To proceed, press Next .

How to verify
Navigation system
is Operational


Figure 62

Step
Verify Navigation System is operational. <i>NOTE: Navigation system connected and Calibration correct icons should both be checked and green. If icons are not checked, verify your Navigation System setup and retry.</i>
Options
To change parameters, press Back .
To proceed, press Next .

6.7.9 Selecting Patient Orientation and Position

How to Select Orientation and Position


Figure 63

Steps
1. Select patient orientation.
2. Select patient position.

NOTE: Display also shows direction of scan.

Options
To change parameters, press Back .
To proceed, press Next .

6.7.10 Selecting 2D Scout Scan View

How to Select the
Scan View

Figure 64

Options

Options	Action
To skip 2D Scout Scan, press Skip 2D Scout Scan .	Perform scout scan steps as described in 6.7.11, 6.7.12, and 6.7.13, but result will be a 3D scan. See 6.7.17.
To change parameters, press Back .	
To proceed, press Next .	

6.7.11 Setting 2D Scout Scan End Position

How to Set the 2D End Position


Figure 65

Step	
Using the left and right arrows a move ring to end position.	
Options	Action
To select the end position	 Press Next d .
To redefine the previous settings	 Press Back s .
To turn on the lasers to aid positioning	 Touch Laser f . <i>NOTE: Touch again to turn Laser off.</i>

Airo produces ionizing radiation. Excessive exposure causes serious injury. Follow the procedures in **2.7 Radiation Safety**.

6.7.12 Setting the 2D Start Position

How to Set the 2D Start Position


Figure 66

Step		
Using the arrows a move ring to start position.		
Options	Action	
To select the start position		Press s Next .
To redefine the end position		Press d Back and reposition the ring.
To turn on the laser to aid positioning		Touch f Laser <i>NOTE: Touch again to turn Laser off.</i>

6.7.13 Verifying 2D Parameters and Scanning

How to Verify Parameters and Scan


Figure 67

Steps
1. Review the displayed scan parameters a .
2. Carry the Pendant Display behind a protective barrier. The CT operator moves behind a protective barrier.
3. Essential personnel move behind a protective barrier. Non-essential personnel exit the room before the scan is conducted and do not re-enter until the ionizing radiation alerts have stopped.

Airo produces ionizing radiation. Excessive exposure causes serious injury. Follow the procedures in 2.7 Radiation Safety.

Options
To scan, press and hold Scan S .
<i>NOTE: The positioning lights show the point of scan on the patient. The scan image is streamed to the display as it is constructed.</i>
To stop the scan at any time, release Scan .
To go back press Back .

To stop the scan at any time, release the Scan button. The data will be saved.

6.7.14 Selecting 3D Scan Region

How to Select a 3D Scan Region


Figure 68

Options

To select a region for the 3D scan, touch and move the boundaries of the highlighted area **a**.

NOTE: Yellow lines indicate width of beam outside the Field of View at isocenter.

6.7.15 Move Ring to Start Position

How to Move Ring to Start Position

Figure 69

Step
Press and hold Scan a to move the ring to the start position.
Options
To reset scan regions, press Back .
To proceed, press Next .

6.7.16 Verifying 3D Parameters and Scanning

How to Verify the 3D Parameters and Scan

Figure 70

Steps	
1.	Review a the displayed scan parameters.
2.	Carry the Pendant Display behind a protective barrier. The CT operator and essential staff move behind protective barriers.
3.	Non-essential personnel exit the room before the scan is conducted and do not re-enter until the ionizing radiation alerts have stopped.
4.	Release Scan S .
5.	To scan, press and hold S Scan . The scan image is displayed as it is reconstructed.

To stop the scan at any time, take your finger off the Scan button. The data will be saved.

6.7.17 Summary of 3D Image Controls

Overview


Figure 71

No.	Component	Function/Usage
a	Thumbnails of main Image View	Changes main image to selection
s	Window/Level	Optimizes image for tissue level selected
d	Image Tools	Manipulate image
f	Next	Send scan to DICOM or USB. See 6.7.22.

6.7.18 Main Image View

How to Change the Main Image View


Figure 72

Step
Touch a thumbnail image. The main image changes to that view.
Options
Axial (transverse)
Sagittal
Coronal

6.7.19 Window/Level

How to Change the Window Level


Figure 73


Figure 74

Steps	
1.	Press a to see menu.
2.	Press S for a selection.
3.	Press a again to hide menu.

Options	Effect
Auto	Automatically chooses a window center and width based on analyzing the intensities in the acquired CT image
Bone	Presets the window center and width values based on the known Hounsfield units of the tissue of interest
Brain	
Lung	
Pelvis	
Custom	Manually adjust the window center and width. - To adjust the value of the window width, drag a finger horizontally across the screen. - To adjust the value of the window center, drag a finger vertically across the screen.

6.7.20 Image Controls

How to Use the Image Controls

Functions	
	To pan
	To zoom in or out
	To rotate
	To scroll up or down
	To set slice thickness of view See 6.7.21 <i>Setting Slice Thickness</i>
	To recenter

Step
Touch a button to activate a function.

Options
To manipulate the image, touch the main image.
To change functions, touch a different button.
To proceed, press Next . See 6.7.22.

6.7.21 Setting Slice Thickness

How to Set the Slice Thickness


Figure 75


Figure 76

Steps	
1.	Press a to see menu.
2.	Press S selection.
3.	Press a again to hide menu.

6.7.22 Sending Scan to DICOM or USB

How to Send the Scan


Figure 77

Steps	Action
To proceed without exporting scan	Press Next . The scan is already stored locally on AIRO.
To send the scan to a DICOM node	1. Select the node. 2. Select Send to DICOM node.
To export to USB	Insert USB drive into USB port on holster. 1. Select Export to USB. AIRO creates a DICOMDIR on the USB and exports the file.
To take another scan or exit scanning	Press Next .
To return to the scan	Press Back .

Options	Action
Anonymize	Send without patient data.
Dosage Report	Include the Dosage Report.
User Friendly Files	For USB export only, adds the .dcm extension to files to allow them to be recognized by extension.

6.7.23 Next Steps

Continue with Same Patient?


Figure 78

Option	Result	See
Yes	Goes to Confirm Ring Rotation screen 	6.7.6
No	Returns to Welcome Screen 	6.6.1

6.8 View an Existing Scan

How to View an Existing Scan

Figure 79


Figure 80

Steps	
1.	From the Welcome screen, select Dicom Viewer. See 6.6.1.
2.	Select a patient name a . <i>NOTE: If the patient name is not on the list, there are no existing scans for that patient on this system.</i>
3.	Select a study s .
4.	Select a series d (a scan identified by time, date, and ID). Press Next .

After viewing the scan:

Options	Action
Send scan to network or USB.	Press Send.
View other scans for this patient or create a new scan for this patient.	Press Home . See 6.7.1.

6.9 Emergency Scanning

6.9.1 Taking Emergency Scans

How to Take Emergency Scans


Figure 81


Figure 82

Steps	See
1. From the Welcome Page, select Emergency Scan a .	
2. Enter patient's weights	
3. Continue scan steps	6.7.2
4. After scanning the patient, use the Edit feature under Existing Scans to replace default name and time stamp with the patient name	6.10

6.10 Editing Patient Data

6.10.1 Editing or Adding Patient Data

How to Edit or Delete Patient Data


Figure 83

Step
Select Dicom Viewer on the Patient Data screen. See 6.7.1
Select Patient name afrom the list or search by patient name or ID.
Select Edit or Delete s in the Patient Name.

Change Patient Information and Save

First Name: DEMO

Last Name: PATIENT89

ID: d714d06b-7a1d-41d1

Gender: ☒ male ☐ female ☐ other

Date of Birth: 1970 Jan 01

YYYY MM DD

Save Cancel

1 2 3 4 5 6 7 8 9 0 - = Back

Tab Q W E R T Y U I O P Enter

Caps Lock A S D F G H K J L

Shift Z X C V B N M . Shift

Space

a

Figure 84

Steps	
1.	Edit patient information as required. <i>NOTE: Emergency scans will have the first name "EMERGENCY" and the last name "SCAN [Date and time]," where date and time are the date and time the scan was taken. AIRO generates and assigns a patient ID.</i>
2.	Select Save to save the patient information to the database or Cancel to ignore changes and return to the Patient Selection screen. a

6.11 After Scanning

6.11.1 Post-Scanning Steps

Where to Find Instructions

Instructions	See
Detach surgical table from column	5.3
Transfer table to column	5.3.2
Put AIRO in transport position	3.3.3
Transport system	3.3.6
Park system	3.3.8

6.12 Service Functions

6.12.1 Service Functions Menu

Select an option from the Service Functions Menu




Figure 85

Options	
	Select Hospital Service a to enter the applications available to Hospital Technicians.
	Select Brainlab Service s to enter navigation configuration settings.

6.12.2 Hospital Service

Hospital Service Main Menu








Figure 86

Options	See
 Select Configure a to export log files, calibrate the Pendant and view the QC pattern.	6.12.3
 Select Add New Phantoms s to add a new test phantom.	6.12.4
 Select Gain Calibration d to perform the daily Gain Calibration.	7.2.3
 Select E-Stop Test f to perform the daily Emergency stop test.	6.12.5


Options	See
 Select Manual Scan g to perform a helical scan selecting a specific kV and mA.	6.12.6
 Select Axial Scan h to perform a scan with no Z translation and a specific kV and mA. Used for dosage testing.	6.12.7
 Select Single Shot Scan j to perform an in place scan without motion. Using a specific kV, mA, angle and duration. Used for dosage testing.	6.12.8
 Select Quality Check Scan k to perform a helical self-check scan of the Gammex 464 phantom	8.3

6.12.3 Option: Configure

How to Select Configuration Options

Figure 87

Options	
	Select Calibrate Touch a to run the Pendant/Touch Screen calibration program. Follow on-screen instructions.
	Select Export logs to USB s to copy log and system configuration file to a USB drive.
	Select QC Pattern d to view the Quality Check image pattern.

6.12.4 Option: Add New Phantom

How to Add a New Phantom


Figure 88

Steps	
1.	Select Add New Phantom from the Hospital Service menu.
2.	Enter a first name for the new Phantom.
3.	Press Next .

6.12.5 Option: E-Stop Test

How to perform a daily Emergency Stop Test


Figure 89

Steps	
1.	Select E-Stop Test from the Hospital Service menu.
2.	Follow on-screen instructions.
3.	Press Next .
4.	Set Start and End point for the scan.

How to perform the E-Stop scan

Figure 90

Steps	
1.	Press and Hold the green Scan button.
2.	Wait until Z-axis translation has started.
3.	Press the red E-stop button.
4.	Observe X-ray beeping tone stops within 0.5 seconds and gantry translation stops immediately.
5.	Press Next to continue or contact service if E-stop test fails.

6.12.6 Option: Manual Scan

How to perform a Manual Scan

Figure 91

Steps	See
1. Select Manual Scan from the Hospital Service menu.	
2. Choose a Phantom from the Phantom Selection list.	
3. Press Next .	
4. Confirm Ring Rotation is correct..	6.7.6
5. Press Next .	
6. Enter desired Manual Scan parameters.(kV, mA, DFOV, Navigation)	
7. Press Next .	
8. Select desired Kernel settings.	
9. Press Next .	
10. Select Patient's Orientation and position.	6.7.9
11. Select and perform 2D Scout scan if desired.	6.7.10
12. Select Start and End points of the scan.	6.7.11
13. Perform the Scan.	

6.12.7 Option: Axial Scan

How to perform an Axial Scan

Figure 92

Steps	See
1. Select Axial Scan from the Hospital Service menu.	
2. Choose a Phantom from the Phantom Selection list.	
3. Press Next .	
4. Confirm Ring Rotation is correct..	6.7.6
5. Press Next .	
6. Enter desired Manual Scan parameters.(kV, mA, DFOV, Navigation)	
7. Press Next .	
8. Select desired Kernel settings.	
9. Press Next .	
10. Select Patient's Orientation and position.	6.7.9
11. Press Next .	
12. Perform the Scan.	

6.12.8 Option: Single Shot Scan

How to perform a
Single Shot Scan


Figure 93

Steps	See
1. Select Single Shot from the Hospital Service menu.	
2. Confirm Ring Rotation is correct..	6.7.6
3. Press Next .	
4. Enter desired Single Shot Scan parameters.(kV, mA, DFOV, Navigation, Angle, Duration)	
5. Press Next .	
6. Perform the Scan.	

6.12.9 BrainLab Service

BrainLab Service Main Menu



Figure 94

NOTE: BrainLab Service functions are to be performed by qualified personnel only.

Options	See
 SAB Navigation Calibration Last Cal: 2/19/2014 4:56:15 PM	Select SAB Navigation Calibration a to calibrate the Curve Navigation system in the Scan Away from Base direction.
 STB Navigation Calibration Last Cal: None	Select STB Navigation Calibration s to calibrate the Curve Navigation system in the Scan Towards Base direction.
 Manual Scan	Select Manual Scan to run the verification portion of the Curve Navigation system calibration process. 6.12.6
 Spiral Phantom Scan	Select Spiral Phantom Scan to run a scan on the Spiral Phantom.

7 MAINTENANCE

7.1 Chapter Overview

7.1.1 Contents

Topics Covered

Section	See
Daily Gain Calibration	7.2
Cleaning, Disinfection, and Sterilization	7.3
Inspection	7.4
Preventative Maintenance	7.5
Temporary Disuse	7.6
Disposal Instructions	7.7

7.2 Daily Gain Calibration

7.2.1 Overview

Calibrate

Calibrations should be performed daily or according to the responsible organization's procedure.

A gain calibration scan consists of a scan of the air. Gain calibrations auto-correct gain and offset settings before diagnostic scans begin.

NOTE: Performance checks are not a substitute for scheduled periodic maintenance performed by an authorized technician.

This section describes the steps for performing gain calibration.

The System requires a verification of the Emergency Stop (E-stop) System through a simulated use test. The System software provides step-by-step instructions on the user interface.

NOTE: In the event of an E-stop failure discontinue use of the System and contact Technical Support.

Definition

Calibration is an image correction procedure to preserve image quality. Calibration corrects for inherent pixel-to-pixel gain variations, non-uniform background or offset signals, and non-linearities in the response of the scanner to the X-ray dose. The non-uniformity of the X-ray beam is corrected to provide a lower noise image across the detector array. [Source: ASTM WK7492 2008]

CT numbers are NOT absolute; misdiagnosis is possible. System and patient variables may effect CT Number accuracy. If you rely solely upon CT numbers without taking the following variables into consideration you could mis-diagnose an image.

The following variable effect CT Number accuracy:

- Partial volume effects of anatomy
- X-ray tube deterioration
- Improperly calibrated system (poorly centered phantom, used wrong phantom, replace current calibration files with extremely old calibration files).
- Beam hardening due to patient anatomy, especially bone.

To reduce CT Number variations:

- Warm-up the X-ray tube whenever the system recommends it; make sure the tube design matches the software configuration parameters.
- Center the patient anatomy of interest in the gantry opening. Select a Field of View that encompasses the patient.
- Acquire comparable images with similar scan and reconstruction choices.

7.2.2 Preparing for Calibration

How to Prepare for Calibration


Steps	See
1. Move AIRO to an area suitable for performing CT scans.	3.3.8
2. Remove the table from the column.	5.3.2

7.2.3 Perform Gain Calibration

Start Up and Enter Scan Mode

Steps	See
1. If necessary, transport and set up for scan.	3
2. Option: Log in.	6.3.2
3. Choose Scan mode.	6.3.3


Select Service Functions from the Welcome Screen

Steps	See
1.  Select Service Functions on the Welcome Screen.	8.3
2.  Select Hospital Service on the Service Functions menu.	8.3
3. Log in.	6.3.2

Select Gain Calibration from the Hospital Service Menu

Figure 95

Step	See
Gain Calibration	
	Select Gain Calibration on the Hospital Service screen 8.3

Confirm Readiness

Figure 96

Steps	
1.	Confirm that the table is removed and the area is safe.
2.	Ensure that the environment is safe for X-ray exposure. <i>NOTE: For instructions to ensure the area meets the safety requirements, see 2.7.</i>
Options	
	To proceed press Next .
	To reset, press Back .

Perform Gain Calibration Scan

Figure 97

Step
Press and hold Green button to Scan .

Pass or Fail Notice After the scan, a notice appears.


Figure 98

Options
AIRO gain calibrations are within specified accuracy. To proceed press Next .
AIRO does not accept gain calibration data. 1.Evaluate the scan area for objects that might interfere with the scan. 2.Remove any objects and re-scan. If AIRO still fails gain calibration after several attempts, mark the system as not for use and contact technical support. See 1.2.1.

7.3 Cleaning, Disinfection, and Sterilization

7.3.1 Cleanable Parts

Warnings and Cautions

The CT system is not waterproof. It is not designed to protect the internal components from the ingress of liquid.

Do not use sprays. Sprays may expose electrical components to moisture which can damage the equipment.

Avoid using phenol-based, corrosive or solvent disinfection agents as they may harm the surface material on the equipment.

Do not use high concentration (>500 ppm) sodium hypochlorite (bleach) solutions as they may cause corrosion damage. Corrosion damage can lead to degradation of structural elements that can result in patient/operator injury.

NOTE: If you are not sure of the properties of a cleaning or disinfectant agent, do not use it.

How to Clean

Clean the system in an area that is acceptable for cleaning and disinfection of the AIRO, and that contains enough room to perform this procedure. To clean the system, you will need the following:

- Super Sani-Cloth Wipes (or equivalent)
- Paper Towels
- Bucket of warm water
- Mild Detergent (dish washing soap or similar)
- Cloth rags or wipes

Steps
1. If not already done, unplug the AIRO from the wall mains.
2. Wear PPE as defined by your facilities bio-safety process.
3. Rotate the Gimbal to the Scan position.
4. Raise the Trumpf Column to its highest position to allow easier wipe down of the bellows.

Use caution when wiping down system as the batteries are still supplying power.

Cleaning

Steps
1. Submerge a cloth rag into the warm water. Ring the rag out so it is not dripping.
2. Wipe down all external surfaces of the Gantry, Gimbal, Base, Trumpf Column, Pendant and Pendant cord, and system power cord. If necessary, use a mild detergent with the water to remove the stains.
3. Once all surfaces have been cleaned, allow to air dry, or wipe dry with a dry rag or paper towel.

Disinfection

Steps	
1.	Use a germicidal disposable wipe (such as PDI Super Sani-Cloth® or equivalent.) to remove dirt or stains. Wipe down all external surfaces of the Gantry, Gimbal, Base, Trumpf Column, Pendant and Pendant cord, and system power cord. The surfaces must stay wet for at least two minutes. Use as many Sani-Cloths as necessary to ensure the surfaces are wet for two minutes. <i>NOTE: The Super Sani-Cloth should not be used to remove dirt and stains. They are for disinfection only.</i>
2.	Allow to air dry.
3.	Wipe the exterior of the bucket and cleanser containers with a Sani-Cloth.
4.	Dispose of all wipes, paper towels, orange sticks and disposable PPE in the bio-hazard container.
5.	Wash hands thoroughly with antibacterial soap.

Interval Clean daily or more frequently as needed, following your organization's sanitation policies.

Tolerance to Water The moisture tolerance of the device is component dependent. Nothing more wet than a wrong out drip-less towel should be used anywhere near the device. The column is rated IPX4. Overall system is IPX0.

7.3.2 Sterilization

No Sterilization	Airo is not intended to be sterilized.

Disclaimers	Brainlab and the system manufacturer bear no responsibility for sterilization of system surfaces. Brainlab and the system manufacturer shall not be held liable for any contagion which may have stemmed from contaminated stains on the system surfaces.

7.4 Inspection

7.4.1 Daily

Wear or Damage The responsible organization and personnel using the machine are responsible for performing daily inspections.

Inspect for:

- Signs of wear or damage

In addition to the above, the system should be evaluated weekly for proper functionality:

- System powers up normally
- Screen display illuminates correctly
- Movement controls operate correctly
- System Air Filter is clean and clear of obstructions.

If upon inspection the system appears damaged or displays functionality inconsistent with the functions described in this manual, contact Technical Support. (See Section 1.2 Contact Data for additional information.)

If the system is used daily, it should be evaluated as described in Weekly inspection, below.

Other Follow the responsible organization's protocols for inspection. See 8 *QUALITY TESTING* for Quality Assurance tests.

7.5 Preventative Maintenance

7.5.1 Preventative Maintenance is Required

Interval Preventative maintenance must be performed once a year to ensure functionality and safety. The owner/responsible organization is responsible for timely maintenance.

Exceeding the one-year period for preventative maintenance is considered improper use.

Authorization Only Brainlab-authorized technicians are qualified to perform maintenance.

Risk of electrical shock. There are no user-serviceable parts. All servicing is to be carried out by Brainlab authorized technicians.

Arrangement

Service Contract	Annual Inspection
Yes	Automatically performed by Brainlab
No	Contact technical support for arrangement

Scope Preventative maintenance includes:

- a complete inspection of the AIRO system
- potential adjustment of dosimetric accuracy
- adjustment or lubrication of normally wearing components

7.6 Temporary Disuse

7.6.1 Temporary Disuse

Definition: Normally, AIRO is left plugged in whenever it is not being moved. Temporary disuse is defined as when AIRO is left unplugged for over 2 days.

Unplugged Over 2 Days Example: The system was accidentally left unplugged over the weekend.

AIRO can be left unplugged for up to one month, but this is not advised.

AIRO batteries require charging, even when the system is unused. Leaving AIRO unplugged causes the batteries to deplete.

Over-discharge can permanently damage the batteries.

How to Maintain Batteries

Leave AIRO plugged in whenever possible.

7.7 Disposal Instructions

7.7.1 Disposal of System and Batteries

How to Dispose of System

Do not dispose of AIRO or its components. Only dispose of electrical and electronic equipment in accordance with statutory regulations. For information regarding the WEEE (Waste Electrical and Electronic Equipment) directive, contact technical support. See *1.2 Contact Data*.

How to Dispose of Batteries

Do not attempt to access or dispose of any batteries in AIRO. The batteries inside AIRO are not accessible by users and potentially hazardous materials. Contact technical support for instructions. See *1.2 Contact Data*.

8 QUALITY TESTING

8.1 Chapter Overview

8.1.1 Contents

Topics Covered

Section	See
Test Object and Software	8.2
Testing Image Quality	8.3
Quality Check Reporting	8.4
Exposure Testing	8.5

8.2 Test Object and Software

8.2.1 Description

Manufacturer and Model Numbers of Phantom and Software

Supplied Quality Assurance Materials	Manufacturer and Model
Test object (phantom)	Gammex 464

Installed Test Image AIRO is delivered from the manufacturer with a phantom scan loaded on the system. This provides a scan for comparison with scans done after delivery and installation.

The phantom description and image are stored with all phantom scans and cannot be deleted.

8.3 Testing Image Quality

8.3.1 Testing Image Quality with the Phantom

Daily Image Quality
Test

Steps	See
1. Scan phantom	8.3.2
2. Evaluate image quality	8.4

8.3.2 Scan a Phantom for Quality Testing

Pre-Scan Checks

Steps
1. AIRO warm-up procedure has been performed.
2. Valid gain calibration scan has been performed.

Start Up and Enter Scan Mode

Steps	See
1. If necessary, transport and set up for scan	3
2. Option: Log in	6.3.2
3. Choose Scan mode	6.3.2

How to perform Quality Check Scan


Figure 99

Steps	See
1. Select Service Functions from the Welcome Menu	6.6
2. Select Hospital Service from the Service Functions Menu.	6.12.1
3. Log in.	6.3.2
4. Select Quality Check Scan from the Hospital Service Menu.	6.12.2
5. Align the Gammex 464 Phantom as shown.	
6. Click Next .	

Perform the Scan

Figure 100

Step
Press and hold the green button a to run a complete scan.

NOTE: AIRO will automatically scan the length of the Phantom, process the data, and run the Quality Checks.

NOTE: Data is automatically and permanently stored on AIRO under Existing Scans with a patient name of Phantom.

8.4 Quality Check Reporting

8.4.1 Examining Daily Image Quality Report

How to view and export the Quality Check Report


Figure 101

Steps	See
1. Use the scroll bar a to view the metrics included in the report.	
2. Insert a USB drive into the available port on the Pendant Holster	3.2
3. Choose Export to USB to copy the report files to a USB drive.	

NOTE: Quality Checks may fail. Ensure the Phantom is correctly aligned against the lasers. If the problem persists contact Service.

NOTE: The image review software on the AIRO scanner does not provide the ability to calculate mean and standard deviation of CT numbers for user-selectable ROIs. However, the software automatically calculates mean and standard deviations of CT numbers for pre-selected ROIs when performing the daily IQ test. This information is presented in the IQ test report and can be viewed using a standard internet browser (Internet Explorer, Chrome, Safari, etc.) via exporting to USB flash drive. For cases where the user needs control over the size and location of the ROI, Mobius recommends importing the scan DICOM data (exported to USB or pushed to a DICOM server) into a 3rd-party application that supports these features.

Quality Check Results Summary Report

Results Summary						
Test Name	Units	Measured	Typical	Minimum	Maximum	Pass/Fail
Alignment Error (X)	mm	0.4	0	-2	2	Pass
Alignment Error (Y)	mm	-0.6	0	-2	2	Pass
Alignment Error (Z)	mm	-0.2	0	-2	2	Pass
Alignment Error (Theta)	deg	-0.1	0	-0.5	0.5	Pass
Alignment Error (Tilt)	deg	-0.03	0	-0.5	0.5	Pass
Alignment Error (Wag)	deg	-0.26	0	-0.75	0.75	Pass
Spatial Linearity (Z-M)	mm/mm	0.9995	1	-0.995	1.005	Pass
Circular Symmetry (C-V)	mm/mm	0.9993	1	-0.997	1.003	Pass
Z Axis Linearity	mm/mm	1.0006	1	-0.995	1.005	Pass
CT Sensitometry (Bone)	HU	939	955	850	970	Pass
CT Sensitometry (Air)	HU	-980	-1000	-1005	-970	Pass
CT Sensitometry (Acrylic)	HU	121	120	110	135	Pass
CT Sensitometry (Poly)	HU	-97	-95	-107	-84	Pass
CT Sensitometry (Water)	HU	2.6	0	-5	5	Pass
SSP FWHM	mm	1.71	1.6	1	1.8	Pass
Tomographic Thickness	mm	1	1	0.95	1.05	Pass
MTF Kernel Sharp (50%)	cycles/cm	2	2	1.6	2.4	Pass
MTF Kernel Sharp (10%)	cycles/cm	3.8	4	3.2	5	Pass
MTF Kernel Sharp (2%)	cycles/cm	6	6	5.5	7	Pass
CT * Accuracy	HU	1.9	0	-5	5	Pass
CT * Noise	HU	4.6	3.5	1.8	5	Pass
CT * Uniformity	HU	2.7	3	0	5	Pass
CT Skew over Z	HU/mm	0.05	0	-0.15	0.15	Pass
CT Standard Deviation over Z	HU	0.3	0.5	0	1	Pass
Contrast to Noise Ratio	HU/HU	2.1	1.2	1	3	Pass

Figure 102

NOTE: The full version of the quality check results can be exported to a USB drive.

•

8.5 Exposure Testing

8.5.1 Responsible Organization

Safety Procedures AIRO meets standards for CT medical equipment cited in 9.12 at point of delivery. After delivery, the responsible organization must comply with the manufacturer's safety and maintenance requirements and follow all standards established by the applicable laws for radiological testing of medical equipment.

The method to measure and verify CTDI on AIRO is provided below.

8.5.2 Measure and Verify CTDI

Equipment Required	1. Polymethyl methacrylate (PMMA) CT dosimetry phantoms - Body phantom 32 cm diameter, head phantom 16 cm diameter, at least 14 cm in length. Examples: Fluke Biomedical 76-415 (body) and 76-414 (head) phantom; Ludlum L-451 (body) and L-441 (head) phantom. 2. Electrometer / pencil ion chamber, 10 cm active detector length, with matched control unit. Examples: Radcal 2026C control unit with 20X6-3CT ion chamber; Radcal Accu-Dose control unit with 10X6-3CT in ion chamber.
Method	1. Place the body phantom on the patient table and position it at the isocenter using the gantry alignment lasers. Ensure that the phantom is centered along the axial, sagittal and coronal planes. 2. Position the pencil ion-chamber in the center hole of the phantom; plug all other holes with acrylic rods (included with the phantom). Set the control unit to read integrated dose, typically in milli-Roentgens. 3. Perform an axial (sequential) scan with no gantry translation, at a typical abdomen technique (example, 120 kV, 200 mAs, N_xT = 32x1.06 = 33.92 mm). Record the exposure displayed on the base unit (mR). Designate the exposure readings as 'Center mR'. 4. Position the pencil ion-chamber at the 12 o'clock position, plugging the center hole with an acrylic rod. Repeat step 3. Designate the exposure readings as 'Entrance mR'. 5. Repeat the same procedure (steps 1-4) with the head phantom, positioned at isocenter. 6. Calculate CTDI 100 as follows, for both body and head phantoms: - CTDI 100, center (mGy) = [Center mR x 100 mm x CF x 0.0088] / [N x T] - CTDI 100, entrance (mGy) = [Entrance mR x 100 mm x CF x 0.0088] / [N x T] where the 100 mm in the numerator is the active length of the pencil ion chamber, - CF is the ion chamber correction factor, if specified by the manufacturer (typically 1.0), 0.0088 is the factor used to convert mR to mGy, - and N x T is the beam thickness in mm (33.92 mm for the Mobius AIRO) 7. Measured CTDI_w (CTDI weighted) as follows, for both body and head phantoms: - CTDI_w (mGy) = [CTDI 100, center x 0.333] + [CTDI 100, entrance x 0.667] Compare the measured value against the calculated value. The calculated value is found by the following formula; - Body: (mA)*(0.0096)*(29.8229) - Head: (mA)*(0.0096)*(39.7604) 8. Compare the measured results against the calculated results and calculate the percentage difference. The percentage difference should fall within 20% and must fall within 30%.

Sources of Variability

1. Variability associated with ion chamber. Generally, properly calibrated ion chambers are expected to perform within an error margin of well under 5%.
2. Variability associated with entrance exposure measurement at 12 o'clock. Since the radiation beam may not turn on and off exactly at the 12 o'clock position during a single spin, some variability may be observed in Entrance mR measurements. If the error in step 7 is found to exceed 20%, it is recommended that the average of measurements at the 3, 6, 9, and 12 o'clock positions is used to calculate Entrance mR, to rule out variability associated with the ion chamber at this location.
3. Operator error. If errors in step 7 exceeds 20% despite ruling out the variability at the 12 o'clock position, precise positioning of the dosimetry phantoms at isocenter, as well as stationary positioning of the gantry with respect to the table during the exposure should be verified.

9 TECHNICAL DESCRIPTION

9.1 Chapter Overview

9.1.1 Contents

Topics Covered

Section	See
Introduction	9.2
Technical Description	9.3
Special Installation Requirements	9.4
Electrical Specifications	9.5
X-ray Specifications	9.6
Performance Characteristics	9.7
Environmental Requirements	9.8
Safety Classifications	9.9
Replacing Fuses and Other Parts	9.10
Network/Data Coupling to Other Equipment	9.11
Regulatory Compliances and Approvals	9.12
RF Communications Equipment	9.13
Symbols	9.14
Labels	9.15
Location of Certification Labels	9.16

9.2 Introduction

9.2.1 Obtaining More Information

**Circuit Diagrams
and Parts Lists**

Brainlab will make available on request circuit diagrams, component part lists, descriptions, calibration instructions or other information that will assist Service Personnel to repair parts that are designated as repairable by Brainlab. See *1.2 Contact Data*.

9.2.2 Warning Regarding Modifications

No User-Modifiable Parts There are no user modifiable parts. Unauthorized modifications may void your warranty or service contract.

Unauthorized modification may result in severe injury or death.

9.2.3 Combining Accessories and Sub-assemblies with AIRO

See Allowed
Equipment List

11 ALLOWED EQUIPMENT

9.3 Technical Description

9.3.1 Brief Description and Operating Principle

AIRO Description

AIRO is a 3rd-generation high resolution, large bore (107cm), 51.2 cm field of view, X-ray Helical Computed Tomography Scanner. AIRO combines a treatment-specific patient support system and an X-ray imaging device that produces high resolution helical images of the body.

It operates by continuous rotation of an X-ray tube and opposing detector elements that irradiate the scanned object from 360 degrees while simultaneously translating with respect to the patient along the patient 'Z' axis following a precise motion profile.

A 'scan' is comprised of a set of 2D axial slice images that are reconstructed from a series acquisition of 2D-projection images. Its design accommodates intra-operative use. The surgical column and patient support are controlled using an integrated keypad or a wireless remote control.

AIRO includes a hand-held control unit for imaging, transport, service, calibration, and remote control used for the integrated table column and associated tabletops.

Unique Features System is mobile: The base incorporates casters, main drive wheel, and an electrical drive system. The gantry is both lightweight and compact. The entire unit is designed to be driven down hallways, on limited inclines, through doorways and onto elevators. The user sees the path of travel through an integrated camera and display. AIRO can be run in any area that has been approved for use by the responsible organization.

NOTE: The Installation Procedure contains the specific requirements for approved areas.

- *Table is detachable and exchangeable:* The column accommodates several models of operating tables and patient supports. The user may change the table and table orientation relative to the base to suit the procedure.
 - *Table is portable:* The table and patient can be moved to the system on a shuttle. The table and patient are then transferred onto the system column. Therefore the operating table can be used in the system.
-

Main Subassemblies AIRO consists of six main subassemblies:

- Base
- Gantry
- Gimbal
- Gimbal Holster and Control Panel
- User Interface Pendant
- Trumpf Column

9.3.2 Significant Characteristics

Significant Characteristics

Attribute	Description
Mobile	AIRO is mobile, and its transport is power-assisted.
Large Imaging Bore	The imaging bore accommodates patients in procedure-specific positions and orientations.
Compact	The bore is aligned with the base for transportation and storage. In transportation mode, AIRO can move through an average sized doorway, be transported by elevator, and traverse a 10 mm (3/8") threshold.
Soft Tissue Imaging	The primary application focus is providing diagnostic-quality imaging of soft-tissue parts of the human anatomy.
Touchscreen User Interface	The CT operator enters information, makes selections, and manipulates the images using a touchscreen interface.
Interchangeable Tables and Patient Supports	AIRO 's column accepts transportable operating tables and other patient supports. The list is available from TRUMPF.
Wireless Remote Control	The surgical column, table, and motion-enabled patient support accessories can be controlled by either a keypad on the column itself or by using a wireless remote control.
Powered by Wall Outlet	AIRO is powered from an ordinary hospital wall 15 amp outlet (AC).
Network/Data Couplings	There are 4 kV isolators on the Ethernet connection.
DICOM Export Capability	Users may export 2D and 3D images in DICOM format to a DICOM service class user (SCU) functioning as a PACS server.
External Memory Devices	Users may also export images and patient data to a USB storage device.
Narrow Gantry Width	Designed to be re-positioned such that the surgeon has clear access to the patient.

9.3.3 How the Device Functions

Modes	AIRO has two main functional modes: <i>Transport Mode</i> and <i>Scan Mode</i> .
Transport Mode	In transport mode the gantry is aligned with the base and the entire system is powered from a battery pack located on the gantry. In this mode AIRO can be driven and positioned by the control pendant attached to one side of the gimbal. Moving AIRO typically requires only a single person even when taking the system up or down a ramp. Once AIRO has been moved to the desired destination (storage area or scanning area) it is powered from line power.
Scan Mode	In scan mode, the base is lowered onto the ground by retracting the casters in the support base. The gantry is rotated 90 degrees perpendicular to the base. In scan mode, all functionality other than transport is available to the user, such as patient loading, system positioning, image acquisition, and image processing. The user interface pendant is used to view and manipulate stored images from the patient database. Images can be archived to the USB digital storage media or pushed to a PACS system through DICOM protocol.
Electrical and Software Architecture	AIRO software converts user commands to actions that AIRO processes to capture, display or transfer images to external devices. When the system is not scanning, the gantry program and the gimbal program exchange data and commands. The software controls the gantry and gimbal hardware using firmware which controls digital and analog control and data signals.

9.3.4 Physical Characteristics

System

Category	Feature	Specification
Physical dimensions	Transport Position:	0.60 x 2.28 x 1.98 m (23 x 90 x 78")
	Scan Position:	1.93 x 3.30 x 1.90 m (76 x 130 x 75")
	Weight:	Without table: 975 kg (2150 lb) With table attached: Varies depending on model and accessories.
	Ring: Exterior diameter:	168 cm (66")
	Bore diameter:	107cm (42")
	Image Field of View:	51.2 cm diameter (20.2")
	Ring Width	38 cm (15")
Transport	Electrical drive system Wheel drive Integrated controls Front-view camera Viewing screen	
Transport motion controls	Speed:	1.6 kph (1 mph)
	Braking:	Typical: <15 cm (6")
Ring (gantry) motion (ranges)	Linear ring (translation – z):	1 m (40")
	Rotation	360° in 1.92 secs.
Applied parts	Operating table or other patient support	
Maximum system load	182 kg (400 lbs)	The maximum load for the patient support may be lower than 182 kg. Do not exceed the maximum load for any applied part.

9.3.5 Diagrams

Ring and Base:
Scan Mode


Figure 103

Ring and Base:
Transport Mode


Figure 104

Scan Position
Footprint


Figure 105

Transport Position
Footprint


Figure 106

9.3.6 Pendant Specifications

AIRO Pendant Specifications


Figure 107


Figure 108

Figure 109

Category	Specification
Physical Dimensions	26 cm (10") x 4 cm (1.6") x 17 cm (7") Weight 1.4 kg (3 lbs)
Touchscreen	4-wire resistive touch display 17 cm (7") x 17 cm (7")
Storage	None
Input Voltage	24 Vdc
Power Consumption	5 W

9.3.7 Operating Table Specifications

AIRO Integrated Operating Table

These are the dimensions of a representative model. Dimensions are given in inches. Actual dimensions vary by the model of operating table being used.


Figure 110


Figure 111

Category	Specification
Patient Weight Limit	Depends on model
Vertical Movement Range	Depends on model
Horizontal Movement Range	Depends on model

9.3.8 Column Specifications

AIRO Integrated Surgical Column


Figure 112

Category	Specification
Height	490 mm to 1040 mm (19" to 41.0")
Weight	150 kg (330 lb)
Lifting mechanism	550 mm (21.7")
Incline of transverse axis (Trendelenburg)	$\pm 45^\circ$
Tilt	$\pm 30^\circ$
Column rotation	0 - 360°
Class of Protection (electrical)	I
Level of protection against ingress	IPX4. Overall system is IPX0.
Mode of Operation	Continuous running with intermittent loading 2 min ON, 8 min OFF

9.3.9 Wireless Remote Control for Column

Wireless Remote Control Specifications


Figure 113

Category	Specification
Dimensions including clip (L x W x H)	(230 x 70 x 55) mm (9.1 x 2.8 x 2.2)"
Weight	approx. 325 g (11.5 oz)
Display	240 x 320 pixels, color
Power supply	Integrated 3.6 V lithium ion battery (rechargeable using charger provided)
IR channel	36 kHz
Wireless transmission frequency range	ISM band, 2.4 GHz
Level of protection against water penetration	IPX4. Overall system is IPX0.

9.3.10 Charger

Wireless Remote Control Charger Specifications


Figure 114

Category	Specification
Dimensions	approx. 175 x 115 x 80 mm (6.9 x 4.5 x 3.2")
Weight	approx. 550 g (19.4 oz)
Operating voltage	100 - 230 V
Frequency	50/60 Hz
Maximum power consumption	110 mA
Charging voltage	5V DC
Charging current	maximum 1A
Class of protection	II
Level of protection against water penetration	IPX4. Overall system is IPX0.
Level of protection when using flammable mixtures in the area	No AP safety
Mode of operation	Continuous load

9.4 Special Installation Requirements

9.4.1 Overview

No User Installation Only qualified Brainlab personnel or their authorized agents may install AIRO. Detailed instructions and requirements are in the Installation Procedure.

Approved Areas Only The responsible organization is to designate areas for use that meet all characteristics specified in the AIRO Pre-Installation Guide Part Number MI-42-0002.

AIRO is to be used only in areas approved for use by the responsible organization. Any use outside of areas that are authorized by the responsible organization is considered improper use.

9.5 Electrical Specifications

9.5.1 Specifications

Electrical Specifications

Category	Specification
Power Source	Internal Suitable for connection to Supply Mains only
Other Power Sources	None
Rated Supply Voltage	100 - 240 VAC
Rated Frequency Range	50/60 Hz
Power Requirements	1500 VA
Phase	Single
Least Favorable Input Power Conditions	100 VAC, 50 Hz
X-ray tube voltage x data values	• Maximum current @ nominal voltage (120 kV): 250 mA • Maximum voltage @ maximum current (250 mA): 125 kV • Voltage, current @ maximum power: 125 kV, 250 mA • Nominal electric power @ 120 kV: 30 kW (120 kV, 250 mA, 20s)
Fuses	Current: 15 A, time delay Voltage: 250 VAC Size: 3 AG (1/4" x 1-1/4")
Voltage Statement	Not applicable. The X-ray system does not operate on line power.
Current Statement	Not applicable. The X-ray system does not operate on line power.
Technique Factors for Maximum Line Current	Not applicable. The line does not power the X-ray system.

Category	Specification
Earth and Supply Terminal Markings	Marked on the system using symbols listed in <i>9.14 Symbols</i> .
Supply Mains Isolation	External (wall plug)
Apparent Resistance of Supply Mains	$R \leq 0.6$ ohms
Proper Length and Gauge of Supply Cable	North America (115V supply mains) • Length: not to exceed 5.0m (16.5ft) • Gauge: 14AWG International (230V supply mains) • Length: not to exceed 5.0m (16.5ft) • Gauge: 1.0mm²

9.6 X-ray Specifications

9.6.1 Generator (High Voltage Power Supply)

Generator (HVPS)

Category	Specification
Type	High Frequency
Location	Gantry
Method of Cooling	Air cooled
Power Rating [kW]	32
Peak power	32 kW
Duty Cycle	2 full energy scans per hour (120 kV x 250 mA x 20 sec, twice per hour)
Input power	Lithium Iron Phosphate LiFePO ₄ battery system

Parameters	Range
Rating	125 kV, 250 mA, 20 s 125 kV, 125 mA, 40 s
Output voltage	80 - 125 kV
Output current	5 - 250 mA
Nominal electric power (1s)	32 kW
Maximum average electric power	521 W
Minimum battery charge for proper operation	280 Vdc
Maximum Operating Points	125 kV, 250 mA, 20 s 125 kV, 125 mA, 40 s 30 minute cooling time between maximum scans
Maximum Deviation from Indications	Tube Voltage: +/-1% Tube current: +/- 2% Tomogram scan times: User-defined Topogram scan times: User-defined Spiral scan times: User-defined

9.6.2 X-ray Source Assembly

X-ray Tube

Category	Specification
Tube type	Rotating Anode
Design	Anode end grounded with metal housing
Beam shape	Cone
Target	7" rhenium-tungsten-molybdenum
Tube voltage and accuracy	80, 100, 120 kV Accuracy: +/- 2%.
Tube current and accuracy	10 - 250 mA Accuracy: +/- 2%
Focal spot size	1.0 mm x 1.0 mm
Total filtration at max kV (inherent + beam shaping filter) on central axis (mm Al equivalent @ 120 kV)	6.8
Anode heat capacity (MHU)	2.3
Maximum anode cooling rate (kHU/min)	2874
Method of tube cooling	50% water, 50% propylene glycol to air

9.6.3 Battery System

Battery System

Category	Specification
Battery type	LiFePO ₄
Nominal terminal voltage	358 VDC
Total battery system energy storage	36080 Wh
Battery Cell Rating Certificate	UL1642 and IEC 62133
Battery Pack Rating Certificate	UN/DOT 38.3

9.7 Performance Characteristics

9.7.1 General Performance Characteristics

Scanning Capabilities

Category	Specification
Rotation time	1.92 sec
Scan range	1 m
Scan time	42 sec for 1 m scan range
Helical scan pitch	1.415
Slice thickness	32 x 1.0 mm
Image reconstruction time	24 images per sec

Image Quality

Category	Specification
Reconstruction matrix	512 x 512

Scanned Projection Radiography

Category	Specification
Maximum SPR length	1 m
SPR dimensions	512 x 1000 mm
Angular position available for SPR	0° and 90°

Connectivity

Gigabit Ethernet
DICOM 3.1

9.7.2 Range, Accuracy, and Precision

Range and Accuracy

Parameter	Value
Range	60 to 125 kVp
Accuracy	+/- 10% over the range

Graph of Accuracy of Radiation Output

Machine Settings		
mA	FS	SDD
25	L	29"


Figure 115

Precision

Parameter	Value
Pass/Fail Criteria	10%
Highest Percent Error	4.7% @ 60 kVp, 25 mA
Lowest Percent Error	1.2% @ 90 kVp, 25 mA

9.8 Environmental Requirements

Environmental Requirements

To avoid damaging the system in transport, storage, or unpacking, observe the marked restrictions on the equipment and any protective packaging.

	Specification	
Category	Operating and Hospital Storage (see note 1)	Transport and Warehouse Storage (see note 2)
Temperature	10°C (50°F) to 25°C (77°F)	0°C (32°F) to 50°C (122°F)
Humidity	20% to 70% non-condensing	5% to 95% condensing
Pressure	70 kPa to 110 kPa	70 kPa to 110 kPa
Local Restrictions	- Do not store near chemical products or gas. - Do not expose to direct UV light. - Protect from exposure to moisture, dust, and salinity.	
Altitude	≤ 3000 m (9840')	Not applicable
Adaptation Time	Minimum of twelve hours after extreme storage conditions	Not applicable
Humidity Preconditioning	Minimum of twelve hours after extreme transport or storage conditions	Not applicable

NOTE:

1. AIRO is installed per the installation procedure.
2. AIRO remains in original transportation packaging.

**Least Favorable
Working Conditions**

The least favorable working conditions are:

Parameter	Least Favorable Condition
Temperature	At or above the high limit of the environmental temperature range: 25°C (77°F). Excessive temperature may degrade image quality.
Humidity	At or above the high limit of the environmental humidity range: 70% non-condensing.
Atmospheric Pressure	At or below the lowest limit of the atmospheric pressure: 70 kPa.

**Least Favorable
Operational
Conditions**

The least favorable working conditions are:

Parameter	Condition
Scans per hour	2
Scan length	0.5 m
Tube loading	120 kV @ 250 mA

9.9 Safety Classifications

List of Safety Classifications

Category	Specification
Type of protection from electrical shock	Class 1
Degree of protection from electrical shock	Table top: B (Applied Part) Rest of system: None
Degree of protection from ingress of water and particulate matter	IPX0
Use in prescence of flammable mixtures	Not suitable for use in an Oxygen Rich Environment of the presence of Flammable Anesthetics
Mode of operation	Non continuous

AIRO is not intended for use within an Oxygen Rich Environment. Risk of Fire.

Do not use AIRO with Flammable Anesthetics. Risk of Fire.

9.10 Replacing Fuses and Other Parts

9.10.1 Replacing Parts

Fuses

Do not attempt to replace parts or otherwise install or maintain the equipment. Only qualified technical service personnel may replace the batteries. Unauthorized service can result in severe injury or death and/or equipment damage.

Fuses Required for AIRO Supply Mains:	
Current	15 A, time delay
Voltage	250 V
Size	3 AG (1/4" x 1-1/4")

The power supply cord is replaceable by authorized technical service personnel only.

Contact technical service for correct replacement of interchangeable or detachable parts that are specified as replaceable (e.g., remote controls).

9.11 Network/Data Coupling to Other Equipment

9.11.1 Overview

Network and Data Coupling	AIRO is designed to be connected to other equipment by network and data couplings.
Characteristics	The necessary characteristics of the network/data coupling for AIRO are Gigabit Ethernet RJ-45.
Hazards	The hazardous situations resulting from a failure of the network/data coupling to provide the specified characteristics are the inability to transfer information or retrieve worklist or patient information.
Instructions	The Responsible Organization is hereby instructed that: • Connecting AIRO to a network/data coupling that includes other equipment could result in previously unidentified risks to patients, operators or third parties. • The Responsible Organization should identify, analyze, evaluate and control these risks. • Subsequent changes to the network/data coupling could introduce new risks and require additional analysis. • Changes to the network/data coupling include: - Changes in the network/data coupling configuration. - Connection of additional items to the network/data coupling. - Disconnecting items from the network/data coupling. - Updating the equipment connected to the network/data coupling. - Upgrading the equipment connected to the network/data coupling.

9.12 Regulatory Compliances and Approvals

9.12.1 Summary of Standards Compliance

List and Description

AIRO was tested to IEC 60601 requirements for ME systems and complies with the following regulatory standards:

Performance Standards	Description
IEC 60601-1:2005	Medical electrical equipment Part 1: General requirements for basic safety and essential performance.
IEC 60601-1-2:2007	Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic Compatibility - Requirements and Tests
IEC 60601-1-3:2008	Medical electrical equipment – Part 1: General Requirements for Safety 3. Collateral Standard: General Requirements for Basic Safety and Essential Performance for Radiation Protection in Diagnostic X-Ray Equipment
IEC 62366:2007	Issued: 2007/10 Ed:1, Medical devices – Application of usability engineering to medical devices
AIRO CT Scanner MobiCT-32 IEC 60601-2-44:2009	Medical electrical equipment – Part 2-44: Particular requirements for the basic safety and essential performance of X-ray equipment for computed tomography
IEC 60825-1:2007	Safety of laser products
FDA 21 CFR 1020.30	Radiological Health – Part 1020: Performance standards for ionizing radiation emitting products – Sec. 30: Diagnostic X-ray systems and their major components
FDA 21 CFR 1020.33	Radiological Health – Part 1020: Performance standards for ionizing radiation emitting products – Sec. 33: Computed tomography (CT) equipment
FDA 21 CFR 1040	Radiological Health Part 1040: Section 10: Performance standards for light-emitting products: Laser products

9.12.2 General Electromagnetic Immunity

Electromagnetic Environment	AIRO is intended for use in the electromagnetic environment specified in the following sections. The owner of the equipment is responsible ensuring that the systems are used in such an environment.
Electromagnetic Immunity Declaration	The tables in the following sections provide guidance according to the manufacturer's electromagnetic immunity declaration.

9.12.3 Electromagnetic Immunity

IEC 61000-4-2,
IEC 61000-4-4,
IEC 61000-4-5,
IEC 61000-4-8,
IEC 61000-4-11
IEC 61000-4-6
IEC 61000-4-3

Immunity Test	IEC 60601 Test Level and Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6 kV contact +/- 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+/- 1 kV differential modes +/- 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$< 5\% U_T$ (>95 % dip in U_T) for 0.5 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $< 5\% U_T$ (>95% dip in U_T) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the AIRO user requires continued operation during power mains interruptions, it is recommended that AIRO be powered from an uninterruptible power supply or a battery.
NOTE: U_T is the AC mains voltage prior to application of the test level.		

Electromagnetic Immunity Testing

AIRO is intended for use in the electromagnetic environment specified below.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Conducted RF IEC 61000-4-6	3 V/rms 150 KHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the AIRO devices including cables, than the recommended separation distance calculated from the equation appropriate to the frequency of the transmitter. Recommend separation distance: $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters as determined by an electromagnetic site survey¹, should be less than the compliance level in each frequency range.² Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

1. Field strengths from fixed transmitters, such as base stations for radio (cellular/ cordless) telephones and land mobile radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the AIRO devices are used exceeds the applicable RF compliance level above, the AIRO device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the AIRO device.

2. Over the frequency range 150 kHz to 80MHz, field strengths should be less than 3 V/m.

9.13 RF Communications Equipment

9.13.1 Overview:

Electromagnetic Environment

Portable and mobile RF communications equipment can affect the systems.

AIRO is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled.

Users can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and **AIRO** as recommended below, according to the maximum output power of the communications equipment.

Interference

NOTE: AIRO should not make any unintended movements due to electromagnetic interference in the presence of high frequency surgical devices, cardiac defibrillators, and their monitors. Interference in movement when using high frequency surgical devices simultaneously with the wireless remote control is possible but uncommon, since no remotely controlled movements are necessary at that time.

Separation Distances

Recommended separation distances between portable and mobile RF communications equipment and **AIRO**:

Output Power of Transmitter (W)	Separation Distance According to Frequency of Transmitter (m)		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. <i>NOTE: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.</i> <i>NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</i>			

9.14 Symbols

9.14.1 Symbols Used

Symbols and Explanations



Symbol	Explanation
	Type B Applied Part according to IEC60601-1
	Warning. Consult the accompanying documents for important safety-related information such as warnings and precautions that cannot, for a variety of reasons, be presented on the device itself.
	Potential equalization point
	Ethernet port for communication with the Image Guidance System software (optional software package from Brainlab).
	Ethernet port for communication with the hospital network.
	Do not reuse
	Non-Sterile
	Do not use if packaging is damaged
	Keep away from sunlight
	Keep dry
	Storage conditions for temperature: The specified temperature range is shown on each label
	Storage conditions for relative humidity non-condensing: The specified humidity range is shown on each label.
	Storage conditions for air pressure: The specified air pressure range is shown on each label

Symbol	Explanation
	Quantity of products in packaging
	Batch number
	Serial number
	Catalog/model number
	Date of manufacture
	Manufacturer
	Authorized representative in the European Community
IPXY	X = Protection against ingress of solid objects Y = Protection against ingress of liquid
	On/off switch
	Standby switch to bring the device into standby mode.
	Emergency Stop
	Follow the instructions for use
	Functional Earth Terminal
	ETL Field Evaluated Mark
	This system complies with the requirements of all relevant standards and directives in effect within the European Community.


9.15 Labels

9.15.1 Labels Used

Labels and Explanations

Label	Explanation
 100-240 VAC 50/60 Hz P.max. 1500VA Mobius Imaging, LLC 323 West Main St. Ayer, MA 01432 USA CE 0123 This equipment is in compliance with the applicable standards of CFR, Title 21, Subchapter J as of the date of manufacture MI-43-0001 Intertek 4009639 04-43-0035 21250672 60601-1 60601-1-3 60601-1-6 62366 60601-2-44 (57), 72 60601-1	System Name Plate/Label: Indicates manufacturer and input power requirements. This equipment is certified to be in compliance with the applicable standards of CFR, Title 21, subchapter J as of the date of manufacture. Location: On the base to the right of the sockets for power and network cables.
 7KLV GHYLFH LV QRW IRU XVH RQ SHGLDWULF SDWLHQW ZKLFK LQFOXGHV SDWLHQW OHVV WKDQ DSSURLPDWHO 50 NJ (110 OEV) LQ ZHLJKW DQG 150 FP (59 LQ) LQ KHLJKW.	Pediatric Caution Label Location: On System gantry above the user interface and holster.
 TRUMPF Medizin Systeme Operating table TruSystem 7500 This device contains RF Module FCC ID: TYOJN5121M0 This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: 1) this device may not cause harmful interference, and 2) this device must accept any interference received, including interference that may cause undesired operation. Software 5.4.1.2 Operating table column TS7500 MOBIUS REF: 1704495 101698132 2012-08	Surgical Column Label
WARNING This x-ray unit may be dangerous to patient and operator unless safe exposure factors, operating instructions and maintenance schedules are observed. MI-43-0004	X-ray Unit Warning

Label	Explanation
	X-ray Tube Label: This label certifies that the X-ray tube meets all applicable federal and international standards for X-ray equipment. Location: On the bottom of the system gantry.
LASER LIGHT. DO NOT STARE INTO BEAM CLASS 2 LASER PRODUCT OUTPUT: 3.5mW, 635nm IEC 60825-1:2007 Complies with 21 CFR 1040.10 and 1040.11, except for deviations pursuant to laser notice No. 50, dated June 24, 2007 M-43-0003	Laser Warning Logotype Label in compliance with IEC 60825-1:2007 The laser product complies with 21CFR 1040.10 and 1040.11, except for deviations pursuant to Laser notice No. 50, dated June 24, 2007. Location: On the base to the left of the sockets for power and network cables.
	Input panel: A - Ethernet port for communication with the Image Guidance System software (optional software package from Brainlab). B - Ethernet port for communication with the hospital network. C - Potential equalization point
	WEEE Label Trash Can
	Ionizing radiation label
	Protective Earth Terminal
	Functional Earth Terminal
	Supply Conductor Terminal

Label	Explanation
 Transporting the AIRO System should only be attempted with the gantry in the transport position. The AIRO System is designed to WUDQVSRUW RQDQ LQFOLQH RI 5C PD[LPLXP. 3OHDVH refer to the AIRO User's Guide for instructions. MI-43-0002	Overbalancing Warning Label
	High Voltage

9.16 Location of Certification Labels

9.16.1 Overview

System Certification Labels


Figure 116

The certification labels **a** are attached to the ring as shown.

Column Certification Labels


Figure 117

The certification labels are attached to the column as shown.

10 CT USER INFORMATION

10.1 Chapter Overview

10.1.1 Contents

Topics Covered

Section	See
No Quick-change X-ray Tubes	10.2
Dosage Factors and Image Quality	10.3
Deterministic Effects	10.4

10.2 No Quick-change X-ray Tubes

10.2.1 No Quick-change X-ray Tubes

Quick-change X-ray Tubes

Airo does not permit the use of quick-change X-ray tubes.

10.3 Dosage Factors and Image Quality

10.3.1 Conditions of Operation

About this Section	This section contains the information for CT users required by FDA 21CFR-1020.33.
kV	The generator is capable of delivering kilo voltages from 60 to 125 kV. The typical kV for both body and head scans is 120 kV. <i>NOTE: Prior to the initiation of a scan, the Pendant Display shows the CTDIvol expressed in units of mGy and DLP in units of mGy.cm, for head and for body standards. It also displays the phantom type and diameter on which CTDIvol values are based. These values are recorded in the scan data.</i>
Scan Time Per Gantry Rotation	In clinical use, the gantry is made to rotate at 30.6 rpm, resulting in a radiation beam-on time per spin of 1.92 sec.
mA	The X-ray tube is capable of current output from 10 to 250 mA, resulting in tube current – time product of 19 to 480 mAs. The typical CT condition of operation for an average size adult is expected to be 192 mAs.
Slice Thickness	$N = 32$ data channels in the z direction (number of tomographic sections produced in a single axial scan). $T = 1.06$ mm acquisition slice thickness (nominal tomographic section thickness) $N \times T = 33.92$ mm (commonly called Beam Collimation or Detector Configuration) Table incrementation, I: At present, this is set at 48 mm per tube rotation (25 mm/sec). Future revisions may provide additional speeds.
Pitch Factor	$N \times T = 48 / 33.92 = 1.415$. Future versions of the system may be provided with additional pitch values.

10.3.2 Dose Information

The following dose information is obtained by using Ludlum L-007N PMMA CT dosimetry phantoms alone on the support platform. Separate dose information is provided for head and body applications.

(i)

The dose information is provided at the following locations:

- a. Along the axis of rotation of the phantom
- b. Along a line parallel to the axis of rotation and 1.0 centimeter interior to the surface of the phantom, with the phantom positioned so that the CTDI is the maximum obtainable at this depth.

The CT conditions of operation are the typical values suggested for CT of the head or body.

CTDI Measured Around Axis of Rotation

CTDI₁₀₀ (center):

Phantom	kV	mA	sec	mAs	Position Center, mGy
Body	120	100	1.92	192	12.54
Head	120	100	1.92	192	33.69

CTDI measured 1 cm below the surface

CTDI₁₀₀ at 1 cm interior to phantom surface, at the location where the exposure is maximum obtainable at this depth. This is the entrance location directly above the center of the phantom, and designated as the 12 o'clock position (if the phantom is viewed from the front of the scanner; i.e. the side of the gantry where the patient table is attached).

Phantom	kV	mA	sec	mAs	Entrance, mGy	Positions
Body	120	100	1.92	192	46.33	12 o'clock
Head	120	100	1.92	192	45.58	12 o'clock

CTDI Measured at Other Angles

CTDI₁₀₀ at positions 90°, 180° and 270° from the 12 o'clock position, designated here as the 3 o'clock, 6 o'clock and 9 o'clock positions, respectively.

Phantom	kV	mA	sec	mAs	Entrance, mGy	Position
Body	120	100	1.92	192	37.48	3 o'clock
	120	100	1.92	192	29.43	6 o'clock
	120	100	1.92	192	33.42	9 o'clock
Head	120	100	1.92	192	41.74	3 o'clock
	120	100	1.92	192	35.47	6 o'clock
	120	100	1.92	192	38.85	9 o'clock

Average CTDI

CTDI₁₀₀ as the average of the 4 entrance locations (CTDI₁₀₀ peripheral):

Phantom	kV	mA	sec	mAs	Peripheral mGy
Body	120	100	1.92	192	36.66
Head	120	100	1.92	192	40.41

**CTDI Per Condition
of Operation (ii)**

CTDI₁₀₀ (center) for different selectable CT conditions of operation. The value at 120kV, 192 mAs is treated as 1.00, and other measurements are normalized to this value. Conditions of operation selected are the minimum, midrange and maximum values of voltage (80, 100, and 120 kV) and tube current - time product (48, 96, 192 and 240).

Phantom	kV	mA	sec	mAs	mGy	Normalized mGy
Body	80	25	1.92	48	0.85	0.07
	100	25	1.92	48	1.80	0.14
	120	25	1.92	48	3.14	0.25
Head	80	25	1.92	48	2.89	0.09
	100	25	1.92	48	5.47	0.16
	120	25	1.92	48	8.42	0.25
Body	80	50	1.92	96	1.71	0.14
	100	50	1.92	96	3.60	0.29
	120	50	1.92	96	6.27	0.50
Head	80	50	1.92	96	5.79	0.17
	100	50	1.92	96	10.93	0.32
	120	50	1.92	96	16.85	0.50
Body	80	100	1.92	192	3.41	0.27
	100	100	1.92	192	7.21	0.57
	120	100	1.92	192	12.54	1.00
Head	80	100	1.92	192	11.58	0.34
	100	100	1.92	192	21.87	0.65
	120	100	1.92	192	33.69	1.00
Body	80	125	1.92	240	4.27	0.34
	100	125	1.92	240	9.01	0.72
	120	125	1.92	240	15.68	1.25
Head	80	125	1.92	240	14.47	0.43
	100	125	1.92	240	27.33	0.81
	120	125	1.92	240	42.12	1.25

CTDI Per Peak Tube Potential CTDI₁₀₀ at the maximum of the 4 entrance locations for different selectable CT conditions of operation. The value at 120 kV, 192 mAs is treated as 1.00, and other measurements are normalized to this value. Conditions of operation selected are the minimum, midrange and maximum values of voltage (80, 100 and 120 kV) and tube current – time product (48, 96, 192 and 240 mAs).

Phantom	kV	mA	sec	mAs	mGy	Normalized mGy
Body	80	25	1.92	48	4.12	0.09
	100	25	1.92	48	7.12	0.15
	120	25	1.92	48	11.58	0.25
Head	80	25	1.92	48	4.42	0.10
	100	25	1.92	48	7.66	0.17
	120	25	1.92	48	11.39	0.25
Body	80	50	1.92	96	8.24	0.18
	100	50	1.92	96	14.24	0.31
	120	50	1.92	96	23.16	0.50
Head	80	50	1.92	96	8.84	0.19
	100	50	1.92	96	15.33	0.34
	120	50	1.92	96	22.79	0.50
Body	80	100	1.92	192	16.48	0.36
	100	100	1.92	192	28.48	0.61
	120	100	1.92	192	46.33	1.00
Head	80	100	1.92	192	17.68	0.39
	100	100	1.92	192	30.65	0.67
	120	100	1.92	192	45.58	1.00
Body	80	125	1.92	240	20.60	0.44
	100	125	1.92	240	35.60	0.77
	120	125	1.92	240	57.91	1.25
Head	80	125	1.92	240	22.10	0.48
	100	125	1.92	240	38.32	0.84
	120	125	1.92	240	56.97	1.25

Dose Profiles Technique Factors

Dose profiles were obtained by exposing Laser Luminescent Dosimeter (LLD) rods obtained from, and analyzed by, Landauer, Inc., Chicago, IL. A plot of the dose profile at isocenter, in free air during one axial scan is presented.

The exposure was made at 120 kV, 25 mA, 2 sec (50 mAs), in the absence of the patient table. The LLD rod, which uses Optically Stimulated Luminescence (OSL) technology, has an active length of 150 mm. Exposure is presented in units of dose (mGy*10⁻²) at intervals of 0.05 mm, resulting in 3000 data points. Thickness of the radiation beam was calculated as 44 mm at full-width half maximum

(FWHM). CTDI₁₀₀ free air was calculated using the integral of the dose profile over the central 100 mm span of the LLD. When extrapolated to 196 mAs, the CTDI₁₀₀ free air obtained is 49.65 mGy.

Dose Profile: CTDI
Free-in-Air
Exposure


Figure 118

Variable	Value				
Free Air at Isocenter	120 kVp	25mA	2 sec	50 mAs	Nominal Beam Width 33.92
Max dose: Mean over position range 57 to 93	936 mGy*10 ⁻²				
Half max dose	468 mGy*10 ⁻²				
Beam Thickness at Full Width Half Maximum (FWHM) position range 57 to 93	44 mm				
Exposure length Integral (central 100 mm)	42965 mGy*10 ⁻² -mm				
Exposure length Integral (central 33.92 mm)	31870 mGy*10 ⁻² -mm				
CTDI (mGy*10 ⁻²)	1267				
CTDI (mGy)	12.667				
CTDI at 192 mAs	48.64 mGy				
Z eff	74.18%				

Geometric
Efficiency in Z
Direction

Exposure at 120 kV, 25 mA, 2 sec (50 mAs)	
Dose integral over 33.92 mm nominal beam width (mGy*10 ⁻²)	31870
Dose integral over 100 mm LLD length (mGy*10 ⁻²)	42965
Efficiency	74.2%

Dose Profile: Head Phantom


Figure 119

Variable	Value				
Head phantom position	Center				
	120 kVp	100 mA	1.92 sec	192 mAs	Beam collimation 33.92
Exposure length Integral	112535.8 mGy*10 ⁻² -mm				
Nominal slice thickness (mm)	33.92				
CTDI (mGy*10 ⁻²)	3318				
CTDI (mGy)	33.18				

**Dose Profile: Body
Phantom**


Figure 120

Variable	Value				
Body phantom position	Center				
	120 kVp	100 mA	1.92 sec	192 mAs	Beam collimation 33.92
Exposure length Integral	42893.3 mGy*10 ⁻² -mm				
Nominal slice thickness	33.92 mm				
CTDI (mGy*10 ⁻²)	1265				
CTDI (mGy)	12.65				

10.3.3 Statement of Maximum Deviations

Maximum Deviations

No value should exceed the allowable deviation of $\pm 10\%$ (suggested). For $CTDI_w$, $CTDI_{vol}$ and DLP, measured values should not differ from displayed values by $\pm 30\%$.

NOTE: The 30% allowable error is suggested because most physicists routinely measure entrance $CTDI_{100}$ at the 12 o'clock position, which can be higher than $CTDI_{100}$ (peripheral) by $>25\%$. Medical Physics testing for ACR accreditation also employs measurement at the 12 o'clock position.

10.3.4 Imaging Performance Information

**Noise, Mean CT#,
and Uniformity –
Head Scan (17cm
Phantom)**

Phantom: Conditions of Operation:	Gammex 464 Uniformity Module.
kVp	120
mA	110
Scan Mode	Helical
Pitch Scan Time	1.41
FOV	21.0cm
Recon Image Pixels	512 x 512
Recon Kernel	Standard
Tomographic Section Length	5.0mm

Typical Values			
Circular ROI Diameters	Mean CT # HU	Noise HU	Uniformity HU
60mm @ ISO	0 +- 5	4	—
15mm @ ISO	0 +- 5	5	—
15mm @ 12 o'clock	0 +- 5	3	0 +- 3
15mm @ 3 o'clock	0 +- 5	3	0 +- 3
15mm @ 6 o'clock	0 +- 5	3	0 +- 3
15mm @ 9 o'clock	0 +- 5	3	0 +- 3

**Noise, Mean CT#,
and Uniformity –
Body Scan (30cm
Phantom)**

Phantom: Conditions of Operation:	Gammex Body Ring and Gammex 464 Uniformity Module.
kVp	120
mA	110
Scan Mode	Helical
Pitch Scan Time	1.41
FOV	51.2cm
Recon Image Pixels	512 x 512
Recon Kernel	Standard
Tomographic Section Length	5.0mm

Typical Values			
Circular ROI Diameters	Mean CT # HU	Noise HU	Uniformity HU
60mm @ ISO	0 +- 7	8	—
15mm @ ISO	0 +- 7	10	---
15mm @ 12 o'clock	0 +- 7	6	0 +- 5
15mm @ 3 o'clock	0 +- 7	6	0 +- 5
15mm @ 6 o'clock	0 +- 7	6	0 +- 5
15mm @ 9 o'clock	0 +- 7	6	0 +- 5

**Modulation
Transfer Function**

Phantom: The CT conditions of Operation	Gammex 464 MTF Module.
kVp	120
mA	110
Scan Mode	Helical
Pitch Scan Time	1.41
FOV	21.0cm
Recon Image Pixels	512 x 512
Recon Kernel	Standard
Tomographic Section Length	1.0mm

Typical Results for Standard Kernel


Figure 121

The deviations in the measure MTF's are within the range of +/- 10%

Tomographic Section Thickness (SSP)

Phantom: the CT conditions of Operation	Gammex464 MTF Module.
kVp	120
mA	110
Scan Mode	Helical
Pitch Scan Time	1.41
FOV	21.0cm
Recon Image Pixels	512 x 512
Recon Kernel	Standard
Tomographic Section Length	1.0mm

Typical 1.5mm FWHM +/- 0.5mm @ 1mm tomographic section length

Maximum Deviations for SSP +/- 0.4mm.

10.3.5 Technical Statements

Tube Housing Assembly

Characteristic	Method of Derivation
Leakage technique factors	< 0.88 mGy/h @ 125 kV/8 mA The observed leakage was below the allowable limit of 1.0 mGy in 1 hour (IEC regulations) and 0.88 mGy in 1 hour (FDA regulations).
Minimum Filtration permanently in the useful beam	4.99 mm Al @ 80 kVp, HVL 4.6 mm Al

Anode Heating and Cooling Curves


Figure 122

Tube Housing Cooling Curve


Figure 123

Tube Rating Charts Maximum allowed tube current (mA) as a function of the following starting heat storage and tube voltages:

Focal spot 1.0 x 1.0 mm
7 degrees, 3 phase, 7200 RPM

Vol. scan time (sec)	Starting Heat storage = 16%			Starting Heat Storage = 33%			Starting Heat Storage = 50%		
1	290 (c)	330	260	290 (c)	330	260	290 (c)	330	260
2	290 (c)	330	260	290 (c)	330	260	290 (c)	330	260
4	290 (c)	330	260	290 (c)	330	260	290 (c)	330	260
6	290 (c)	330	260	290 (c)	330	260	290 (c)	330	260
10	290 (c)	330	260	290 (c)	330	260	290 (c)	330	260
20	290 (c)	330	260	290 (c)	330	260	290 (c)	330	260
30	290 (c)	330	260	290 (c)	330	260	290 (c)	250 (a)	200 (a)
40	290 (c)	330	250	290 (c)	260 (a)	210 (a)	240 (a)	190 (a)	150 (a)
50	290 (c)	270 (a)	210 (a)	260 (c)	210 (a)	170 (a)	190 (a)	150 (a)	120 (a)
60	280 (a)	230 (a)	180 (a)	220 (a)	180 (a)	140 (a)	160 (a)	130 (a)	100 (a)

Key	
(a)	Limited by available heat storage
(b)	Limited by window heating
(c)	Limited by filament emission

X-ray Controls and Generators

Parameters		Range	Accuracy
Generator	Power	31.25 kW	Maximum
	Reproducibility	60 - 125 kVp	+/-10% (CV) of 0.05
	Linearity	25 - 125 mAs	Adjacent: (CV) of 0.05 Non-adjacent: ≤0.10
	Rating and Duty Cycle	125 kV, 250 mA, 20 s	
	Output voltage	80 - 125kV	
	Output current	0 - 250 mA	
	Nominal electric power (1s)	326W	
	Maximum average electric power	500W	
Minimum battery charge for proper operation	280 Vdc		
Maximum Operating Points	125 kV, 250mA, 40 s 20 minute cooling time between maximum scans		
Maximum Deviation from Indications	Tube Voltage: +/- 1% Tube current: +/- 2% Tomogram scan times: User-defined Topogram scan times: User-defined Spiral scan times: User-defined		
Exposure Mode (continuous)	Time		
	mA		
	kVp		
Statement of Reference Loading Conditions	Continuous anode input power: 1000 W Nominal X-ray tube voltage: 120 kVp Continuous current under this condition: 8.3mA.		

10.3.6 Definitions of the Measurement Bases of the Technique Factors

Tube Voltage	The kVp can be measured software assisted during each scan via built-in HV-dividers. The values have been verified at a prototype test by using additional devices like a hyper-pure Ge-detector and an electronic penetrameter.
Tube Current	The mA have been verified at a prototype test using a multimeter connected to the center ground point of the high tension transformer. The tube current can be measured software assisted during each scan.
Scan Time	The scan time is the total exposure time (X-ray on) during a scan. It can be measured software assisted during each scan and is verified at a prototype test measuring the X-ray on-signal with an oscilloscope.

10.3.7 60601 Requirements: CTUSERS

60601 IEC 60601-2-44 requires the following information in addition to requirements already met in this chapter.

Nominal Focal Spot Value

1.0 x 1.0 mm

Effect of Radiation Detector on Image

The radiation detector is used to normalize the detector data. This minimizes artifacts due to variations in intensity of the x-ray source.

Half-Value Layers

Half layer values (HLV) were measured using a Radcal 2026 with its 20X6-6 ion chamber placed at isocenter.

Parameter	Value		
Tube Current	25 mA		
Beam-on Time	2 sec (50 mAs)		
kV	80	100	120
mm Al	mR		
0	489	811	1182
3			838
4	264.4		755
5	231.4	435	
6		391	629
7			
8			522
Half exposure	244.5	405.5	591
Measured H _v L, mm Al	4.6	5.7	6.7
Minimum Required, mm Al	2.9	3.6	4.3

Filter Properties

4.99 mm Al @ 80 kVp, HVL 4.6 mm Al

There are no other user selectable filters on AIRO, including shaped filters for head protocols.

Protection Against Residual Radiation

The gantry-side radiologic exposure is 5.7% and 8% of the corresponding exposures on the axis of rotation, thereby meeting the IEC 60601-11 requirement of <25%

10.4 Deterministic Effects

10.4.1 Radiation Dose Levels

Managing Radiation Dose Levels

There is a possibility in normal use that the patient can be exposed to radiation dose levels where deterministic effects may occur. Operators need to be aware of and manage the possibility of exposure to high radiation.

Deterministic Effect: 1000 mGy CTDI ₁₀₀ (peripheral)	Number of Scans to Reach Effect
Typical/maximum head protocol: 120 kV, 326 effective mAs	10
Typical body protocol: 120 kV, 149 effective mAs	25
Maximum body protocol: 120 kV, 326 effective mAs (e.g., obese patient)	11

10.4.2 Prevailing Radiation Dose

Factors Affecting Prevailing Radiation Dose

AIRO selects a protocol based on the data entered for:

- Patient weight
- Body region being scanned

The operator sets the patient weight, anatomy, and reconstruction slice thickness, all of which affect the radiation dose in normal use. See 6 *SCANNING*.

AIRO calculates and displays the kV and mA for the proposed scan.

The Pendant Display shows the CTDIvol expressed in units of mGy and DLP in units of mGy.cm, for head and for body standards. It also displays the phantom type and diameter on which CTDIvol values are based. These values are recorded in the scan data.

NOTE: These clinical protocols are proposed by the manufacturer, pre-loaded on the equipment, displayed prior to scan, and recorded with each scan. They constitute recommendations to be applied directly so as to allow optimized operation. See the Mobius Imaging LLC Application Guide for Protocols and Principles (document MI-42-0005).

11 ALLOWED EQUIPMENT

11.1 Chapter Overview

11.1.1 Contents

Topics Covered

Section	See
Compatibility with Other Medical Devices and Software	11.2

11.2 Compatibility with Other Medical Devices and Software

11.2.1 Medical Devices

Accessories are Marked

All components and accessories that can be removed in normal use are labeled on **Airo**. Information on the label includes a part number, model number, and the name or trademark of the manufacturer to reference for obtaining replacements.

Other Medical Devices

Using medical device combinations which have not been authorized by the manufacturer may adversely affect safety and/or effectiveness of the devices and endanger safety of patient, user and/or environment. Use of unauthorized devices is improper use.

Brainlab Navigation Systems

Brainlab Navigation Systems. Contact Brainlab for more information for proper connection and operation with the AIRO System.

11.2.2 Medical Software

Authorization

Only authorized representatives may install software on the AIRO system. Do not install or remove any software applications.

11.2.3 Accessories List

List of Allowed Accessories

Accessories include:

- Wireless remote control and charging stand for operating the surgical column.
- Trumpf Table tops . Only Trumpf Radiolucent Table Tops (360 Radiolucent only e.g., no metal rails) should be used with the AIRO and integrated Trumpf Table Column.

NOTE: No metal objects should be in the area of the scan. Mobius is not responsible for the quality of scans when metal interference is introduced.

Table tops marked as “Fully Radiolucent” are not sufficient for CT scanning. Only 360 radiolucent tables and table sections should be used in conjunction with AIRO.

Contact Trumpf to purchase associated Tabletops at: (860) 255-6324

12 TROUBLESHOOTING

12.1 Chapter Overview

12.1.1 Contents

Topics Covered

Section	See
Troubleshooting	12.2
Metal Artifacts in Scans	12.3
Transport Failure	12.4
Table Transfer Errors	12.5
Table Malfunctions	12.6
Malfunctions and Return Instructions	12.7

12.2 Troubleshooting

12.2.1 Problems and Suggested Actions

Possible Problems

Problem	Possible Cause	Steps
Cannot connect to hospital network	Network cable is not connected	Check both ends of cable for proper connection try installing a new cable
Metal artifacts in CT scans	Metal outside scan field of view	• Be sure you are using the proper Trumpf Radiolucent 360 patient supports. See <i>11.2.3 Accessories List</i>. • Re-evaluate patient, surroundings, and accessories for metal
	Metal in patient	Use soft kernel setting
Wireless remote cannot find table	Out of range	Move closer. Confirm remote is charged and that this remote is paired to the system.
Transport failure	The drive mechanism failed.	See Section 12.4.
	The mechanism to raise the base failed, or the system can't transition from scan mode to transport mode.	
AIRO will not setup for scan mode.	Internal operation issue.	Restart AIRO, if problem persists contact service.
AIRO will not power on.	Power source interrupted.	Check to make sure AIRO is plugged in, restart AIRO, If problem persists contact service.
AIRO experiences complete loss of power and needs to be moved.	Power source interrupted.	See Section 12.4.
During transport AIRO becomes hung up on obstacle.	Operator error.	Check under machine and remove obstacle if possible. Move the machine in an opposite direction of travel then when the obstacle was encountered. Hold down the forward or reverse button and push the machine (may take a couple of people) in the same direction as the button press would move the machine if under power. See instructions for complete loss of power.
Gain Calibration attempted multiple times, will not pass.	Object inside ring scan area.	Make sure Ring is clear of all objects, If problem persists contact service.
Pendant touch screen does not respond.	Pendant hardware damage or cable connection issue.	Check for damage to Pendant, restart AIRO, If problem persists contact service.

Problem	Possible Cause	Steps
Pendant buttons do not respond.	Pendant hardware damage	Check to see if the AIRO Pendant buttons stuck in the "on/pressed" position . If damage is evident contact service.
AIRO will not power on, Trumpf column needs to be moved.	Power interruption.	Requires two Users: User 1: Turn AIRO Power key on, then press an hold start button on control panel. Hold for 5 seconds and continue holding. User 2: Turn the power to the Trumpf column and manipulate the table. User 1: Release start button and system will power off.

12.3 Metal Artifacts in Scans

12.3.1 Recognizing Metal Artifacts in CT Scans

How to Recognize Metal Artifacts in Scans

Figure 124

Metal objects in scans create streaks that degrade the image quality, and make the scan difficult to interpret. X-ray beams that pass through or near metal implants are highly attenuated and have a much larger error.

NOTE: Image provided for illustrative purposes only.

How to Avoid Metal Artifacts in Scans

Likely Cause	Steps
Metal outside scan field of view	• Be sure you are using the proper Trumpf Radiolucent 360 patient supports. See <i>11.2.3 Accessories List</i>. • Re-evaluate patient, surroundings, and accessories for metal
Metal in patient	Use soft kernel setting

12.4 Transport Failure

12.4.1 Transport Failure

What to Do if You Cannot Move AIRO

Step
Evaluate the situation. Is AIRO raised onto its wheels, or lowered onto the ground? Is the system unplugged? (You must unplug the AIRO before transport.)

Situation	Likely Cause	Response
- The system is raised on its casters. - You are using the pendant, but the system will not move.	The drive mechanism failed.	See 12.4.2
- The base is on the floor. - You cannot raise the system to move it.	The mechanism to raise the base failed, or the system can't transition from scan mode to transport mode.	See 12.4.3

12.4.2 Base is Lifted but System Will Not Move

How to Move AIRO
When the Drive
Wheel Fails

Situation	Solution
The system is on its casters and you are using the pendant, but the system will not move.	Put the main drive wheel into neutral manually. Push from the normal position.

Never attempt to put the main drive wheel into neutral on an incline. There is a risk of injury and equipment damage from uncontrolled movement.


Figure 125

Figure 126

Steps
1. Lift a latch to release gimbal.
2. Rotate the gimbal out of transport position into scan position. This reveals s the drive mechanism in the base.
3. Depress d the release pin. This puts the main drive wheel into neutral.
4. Rotate the gimbal to transport position and lock.
5. Push AIRO manually.
6. Mark the system as not for use and contact technical support.

12.4.3 Cannot Raise System

How to Move AIRO When the Base Will Not Lift This manual method of moving the system is restricted to short distances and level surfaces only

Situation	Solution
The system is lowered and will not rise onto the casters.	Partially raise the base manually. Push from the normal position.

Never attempt to raise the base on an incline. There is a risk of injury and equipment damage from uncontrolled movement.

Equipment Needed: The emergency kit is located under panel covering the drive mechanism and consists of:

Emergency Kit

- 6 mm allen wrench
- 12 mm ball screw
- Phillips head screw driver (The screw driver is not under the panel as it is used to remove the panel.)

Remove Cover


Figure 127

Steps
1. At the top of the ring, tilt the ring by pushing it firmly. (This step may require two people.)
2. Remove the four screws from the panel covering the drive mechanism.
3. Remove the panel.

Raise the System

Steps	
1.	Insert the screw from the emergency kit into the center holes ball and down. (See page 261.)
2.	Using the allen wrench, tighten the screw until AIRO begins to lift.
3.	Turn until the machine lifts about 1/4" off the ground (approximately four turns).

How to Transport

Steps	
1.	Rotate the gimbal into the locked position.
2.	Lock the gimbal in transport position.
3.	Push system from the usual position to move it.
4.	Label the system as not for use and contact technical support.

12.5 Table Transfer Errors

12.5.1 Transfer Errors

Error: Tabletop Direction Undefined

In exceptional cases, the column sensors might not detect the tabletop position. In this case:


Figure 128

- The remote control displays **Direction Table Top undefined, please select manually**.
- The **a Malfunction** symbol illuminates on the column keypad.
- Both **S d** head position indicators blink at the same time.

How to Manually Adjust the Head Position

Manually adjust the head position by pressing the appropriate key.

Options
To set the head position left, press s .
To set the head position right, press d .

126 Table Malfunctions

12.6.1 Situation: Unable to Remove Table

If the Table Will Not Come Off the Column There are two conditions in which the table cannot be removed from the column using the normal functions. **Airo** has features to override these malfunctions:

Reason for Malfunction	See
Column will not return to zero position (drive failure)	12.6.2
Table will not release (sensor failure)	12.6.3

Exercise caution. Remove patients manually if necessary.

12.6.2 Column Fails to Reach Zero Position

Definition The emergency transfer function allows you to remove a table from the column even if the **Zero Position** function fails.

What Happens When the Zero Position Function Fails


Figure 129

Steps	
1.	Press the Transfer pedal a or Zero Position key s to start removing the table from the column.
2.	Zero position halts before reaching zero position. The remote control displays Initialization of [defective drive] not possible! Press key again to continue.
3.	Press the key or the pedal again. The remaining working drives travel to zero position. • A signal (three high-pitched tones) audibly reports the incorrect position. • The remote control displays Incomplete Initialization of Position. • The Zero Position d indicator on the status bar is orange. • Zero Position f blinks on the column keypad and Malfunction g illuminates.

How to Complete the Transfer

Collisions can cause material damage or patient injury.

Section segments on defective joints may not move up and may collide with the shuttle during transfer. Monitor transfer and prevent collisions by aborting the function.

Steps	
1.	Release the Transfer pedal A or Zero Position key S .
2.	The remote control displays Request emergency transfer, with 0-key or pedal .
3.	Press the Transfer pedal A on the shuttle or the Zero Position key S on the remote control for 5 seconds and then release.
4.	Press and hold down the Transfer pedal A on the shuttle or the Zero Position key S on the remote control again. The remote control displays Emergency transfer active. Attention! Check for collision!
5.	Airo completes the transfer at half speed. Airo emits three high-pitched tones (warning signal) continuously
6.	When the column is in the lowest transfer position and the warning signal stops, the transfer is complete. Move shuttle away.

If Airo Will Not Initiate an Emergency Transfer

If the vertical drive is defective, Emergency transfer is not possible.

Step
Move the patient manually.

12.6.3 Table Will Not Release

If the Table Will Not Release If the table will not release due to a sensor defect, emergency operation must be activated.

NOTE: Activate the operating table emergency mode only when an error is present. Emergency mode is only intended for moving the patient.

*NOTE: During emergency operation, **Airo** does not detect when zero position and shuttle positioning are reached. The user evaluates the status.*

How to Complete the Transfer

Figure 130

Steps	
1.	Press On a and Emergency Operation s on the column keypad simultaneously for approx. 5 seconds. You will hear various acoustic signals. Emergency Operation will continually blink red
2.	Bring the operating table into zero position. If necessary, adjust the functions individually.
3.	Run the shuttle up to the stop on the operating table column (no sensor detection)
4.	Use the Lower function to transfer the operating tabletop to the shuttle.

If Airo Will Not Initiate an Emergency Transfer

If Emergency transfer is not possible:

Step
Move the patient manually.

Next Step:
Contact Customer Service

Complete the current procedure and then contact technical support. Prohibit use of the defective operating table for further procedures.

12.7 Malfunctions and Return Instructions

12.7.1 Malfunctions

System Damage or Failure

Do not continue to use the system if:

- The power cable or plug is damaged (e.g., has been run over) or frayed
 - The system does not operate normally when operating instructions are followed.
-

How to Respond to Damage or Failure

Steps	
1.	Turn off system.
2.	Unplug system from wall outlet.
3.	Contact technical support.
4.	Attach a notice such as "DO NOT USE" to equipment to prevent it from being used inadvertently.

12.7.2 Return Instructions

Reporting Damaged Equipment	Any defective components should be immediately reported to technical support.

Removing Components	Only remove defective components if instructed to do so by technical support.

How to Return Components	Follow the instructions from technical support.

How to Contact Technical Support	See <i>1.2.1 Contacting the Company</i> .

INDEX

Numerics

3D Image Controls	142
60601-6 classifications	25

A

Accessories	
combining	193
Add New Phantom	155
Airo	
general overview	4
medical purpose	6
technical description	190
Annual inspection	178
Approved areas for use	206
Assessing patient risk	27
Audio signals	65
Automatic shutdown	
table and column	80
Axial Scan	156

B

Battery charge	
minimum for operation	209
Battery preservation	
automatic shutdown	80
Beam shape	210
Before Scanning	111
Brainlab support	2
Brake (table)	82
Brakes	51
Bypass 2D scan	135

C

Cables	
power and network	44
Calibrate	166
calibration instructions	191
Cautions	4, 7
CE Label	10
center of gravity	102
Change Patient Information	152
Charging remote control	72
Circuit diagrams	191
Collisions	
avoiding	34
collisions	104
Column	63
automatic shutdown	80
controls	66
diagram	64
lock and unlock two-key function	81
power controls	80

rotate	82
Column keypad	66
diagram	67
visual signals	70
Column specifications	203
Compliance with standards	219
Compliances	218
general electromagnetic immunity	219
Component part lists	191
Components	
overview of	8
Confirm the Ring Rotation	130
Connecting remote control	71
Connectivity	212
Contact Data	2
Contact data	2
Contraindications	5
Control	
place of	7
Controls	42
CPR	18
Current	
tube	210
Current statement	207
Custom configurations (remote control)	99

D

Deviation from indications	
maximum	209
Device description	
general	5
technical	194
Diagram	
column	64
column keypad	67
operating table	202
pendant display	201
remote control	73
sample OR setup	54
scan mode	9
transport mode	8
Diagram of System	109
Diagrams	199
with dimensions	199
DICOM	196
conformance statement	10
DICOM 3.1	212
displayed scan parameters	138
Distributor	2
Documentation	13
available guides	13
intended audience	13
requirements	11
Dosages	
preset	28
Driving the system	41

E

Electrical	
safety warnings and cautions	35
specifications	207
Electromagnetic immunity	220
Emergency Scan	122
Emergency Scans	150
Emergency Stop	
location of button	42
next steps	17
Emergency stop	16, 225
releasing patient after	17
when to press	17
EMI	221
End Position	136
Environmental	214
Environmental requirements	214
E-stop test	114
Ethernet	212
Exam Region	131
Existing Patients	125
Existing Scans	121
External memory devices	196

F

FDA 21 CFR 1020.30	219
FDA 21 CFR 1040.10	219
Field of View	118
Field of view	118
Fire safety	37
Flex	89
Focal spot size	210
Footprint	200
Frequency range	207
Functional earth terminal	227
Fuses	207, 217
specifications for	217

G

Gain Calibration	114
Gantry and base	199
diagram with dimensions	199
General Precautions	102
Gigabit ethernet	212
Graph	213

H

Hand grips	42
Hand position	51
Head position	85
override	85
Helical scan pitch	212
Holster	42
hospital network	125
Hospital Service	155
Hospital Worklist	128
Housing type	210
How to Log In	112
Humidity preconditioning	214

I

IEC	219
IEC 61000	221
Image	
quality	212
reconstruction time	212
Implanted medical devices	28
Improper use	
defined	7
Incorrectly entered scan	30
Indications for use	5
Inspections	
annual	178
weekly	178
Installation requirements	206
Intended use	6
Intended users	6
interrupt a transfer	107

K

Key switch	42
Keypad (remote control)	73, 76

L

Label	
by location	229
functional earth terminal	227
high voltage	228
input panel	227
ionizing radiation	227
laser aperture	226
laser warning	227
location of	229
overbalancing warning	228
protective earth terminal	227
supply conductor terminal	227
system rating plate	226
WEEE	227
X-ray tube	227
X-ray unit warning	226
Laser safety	29
Legal information	10
Light	92
Limited Functions	102
Load capacity	
patient weight	6
table	6
Location	6

M

Main Image View	143
Maintenance	165
Maintenance personnel	
qualifications for	192
Maintenance Scans	114
Malfunction	37
Manual Scan	161
Maximum deviation from indications	209
Maximum operating points	209
Maximum patient weight	6
Mechanical motion safety	32
Medical devices	

body worn	28
Medical purpose	6
Metal objects	120
Mobius Imaging LLC	2
Modifications	192

N

navigation system	125
Network/data couplings	196
New Patient	125
New Scan	121
Non-Essential Objects	119
Notes	4

O

Omit 2D scan	135
On	80
Operating points maximum	209
Operating principle	194
Operating room placing Airo in	54
Operating table specifications	202
options menu	122

P

Parking	59
Patents	10
Patient CPR	18
radiation risk	27
Patient Environment	33
Patient environment defined	33
definition of	33
safety	34
Patient Group	6
Pendant Display	109
Pendant display	42
Pendant holster	42
Performance Characteristics	212
Performance standards IEC	219
Phantom	181
Phase	207
Physical Characteristics	198
Position numbering	95
save current	95
selecting stored	95
stored	96
Position menu	94
Power controls	42
current statement	207
frequency range	207
fuses	207
least favorable input conditions	207
maximum line current	207
on and off	41
rated frequency range	207
rated supply voltage	207
requirements	207

source	207
supply cords	35
supply mains isolation	207
surgical column	80
switch	42
terminal markings	207
tube voltage	207
turn on	41
turn on system	45
voltage statement	207
Powering down	59
Preconditioning	214
Preset dosages	28
Profile of stray radiation	25
protective barriers	25
Protective clothing and devices	25

Q

Qualifications maintenance personnel	177
users	6
Quality Check Scan	156
Quality Check Scan (SAB)	114
Quality Check Scan (STB)	114

R

Radiation patient risk	27
profile of stray	23
protective barriers	25
protective clothing	25
safety	25
Radiation safety	20
Radiological safeguards	207
Range	213
Reconstruction matrix	212
Reconstruction time	212
Reflex	92
Regulatory compliances	218
Remote control	66, 71
charger specifications	205
charging	72
functions specific to	90
instructions for using	73
keypad	73
light	93
position menus	94
slow	91
specifications	204
status bar	73
touchscreen	73
Replacing parts	217
authorization	217
other parts	217
Requirements documentation	11
environmental	214
installation	206
qualifications for operators/users	6
training	11
Restart a Stopped Transfer	103
Restrictions least favorable conditions	215
shipping	214
United States	10
Risk reduction	14

Rotate	
column	82
table	82
Rotation time	212
rotation time	212

S

Safety	
fire	37
hazard notices	16
mechanical motion	32
patient environment	34
transportation	36
ventilation	55
warnings and cautions	14
Scan	
pitch	212
range	212
reconstruction time	212
slice thickness	212
stopping unintended	17
Scan field of view	118
Scan Mode	113
Scan mode	9
Scan time	212
Scanned Projection Radiography (SPR)	212
Scanning	
performance characteristics	212
scatter plot	23
Scout Scan	135
Select a Navigated Scan	132
Send the Scan	147
Separation distances for RF communications equipment . 223, 224	
Service (remote control menu item)	97
Service Functions	121, 154
Single Shot Scan	156
Skip 2D scan	135
Slice Thickness	146
Slice thickness	212
Slow	91
Source assembly specifications	210
Specifications	
electrical	207
Standards compliance	219
Start Position	137
Startup sequence	42
Status bar (remote control)	73, 74
Stopping	51
Storage	
environmental requirements	214
Stored positions (table)	96
Stray radiation	23, 25
Supply mains isolation	208
Supply voltage	207
Support	2
Surgical column	63
rating plate	226
Surgical setup	54
Surgical table	
controls	78
Symbols	224
warnings and cautions in user guide	4
System Lock	123
System rating plate label	226

T

Table	
automatic shutdown	80
automatic stops	89
controls	79
custom configurations	99
rotating	82
Table functions	79
all other	89
custom	99
flex and reflex	92
flex/reflex	89
lift and lower	83
light	93
lock and unlock two-key operation	81
on and off	80
other	89
position	94
rotate and brake	82
select head position	85
service	97
slow	91
tilt	86
Trendelenburg and anti-Trendelenburg	87
zero position	84
Table is Not Locked	104
Task Bar icons	115
Technical description	190
Technical support	2
Technique factors	
maximum line current	207
Terminal markings	208
Test Object and Software	181
Tilt (table)	86
Touchscreen (remote control)	73, 75
Training	11
Transfer (zero) position	84
Transferring Patients	101
Translate the Ring	117
Transport mode	41
diagram of	8
Transportation safety	36
Transporting equipment	41
Trendelenburg and anti-Trendelenburg	87
Tube	
current	210
type	210
voltage	210
Two-key function	
column keypad	81
disabling	81

U

USB port	42
User	3

V

Ventilation	55
Visual signals	
column keypad	70
Voltage	207
statement	207
tube	210

W

Warm-up Scan.....	114
Warnings	4
WEEE	10, 179
Welcome Menu	121
Window Level	144
Wireless remote control.....	73, 196
Working conditions least favorable	215

X


X-ray	
field orientation.....	212
generator assembly specifications.....	211
source assembly specifications.....	210
specifications	209
voltage and data values.....	207

Z

Zero position	84, 103

Notary Public hereby certify that this document is a true copy of a

dated September 24, 2015

Dr. Bernhard Schaub

Notary Public in Munich

registered office in 80331 Munich, Tal 12, Germany
commission expires: October 31, 2032

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Results Summary						
Test Name	Units	Measured	Typical	Minimum	Maximum	Pass/Fail
Alignment Error (X)	mm	0.4	0	-2	2	Pass
Alignment Error (Y)	mm	-0.6	0	-2	2	Pass
Alignment Error (Z)	mm	-0.2	0	-2	2	Pass
Alignment Error (Theta)	deg	-0.1	0	-0.5	0.5	Pass
Alignment Error (Tilt)	deg	-0.03	0	-0.5	0.5	Pass
Alignment Error (Wag)	deg	-0.24	0	-0.75	0.75	Pass
Spatial Linearity (LM)	mm/mm	0.9995	1	-0.995	1.005	Pass
Circular Symmetry (XY)	mm/mm	0.9993	1	-0.997	1.003	Pass
Z Axis Linearity	mm/mm	1.0006	1	-0.995	1.005	Pass
CT Sensitometry (Bone)	HU	939	955	950	970	Pass
CT Sensitometry (Air)	HU	-990	-1000	-1005	-970	Pass
CT Sensitometry (Acrylic)	HU	121	120	110	135	Pass
CT Sensitometry (Poly)	HU	-97	-95	-105	-84	Pass
CT Sensitometry (Water)	HU	2.6	0	-4	4	Pass
SSP FWHM	mm	1.71	1.6	1	1.8	Pass
Tomographic Thickness	mm	1	1	0.85	1.05	Pass
MTF Kernel Sharp (50%)	cyc/cm	2	2	1.6	2.4	Pass
MTF Kernel Sharp (10%)	cyc/cm	3.8	4	3.2	5	Pass
MTF Kernel Sharp (2%)	cyc/cm	6	6	5.5	7	Pass
CT # Accuracy	HU	1.9	0	-5	5	Pass
CT # Noise	HU	4.6	3.5	1.8	5	Pass
CT # Uniformity	HU	2.7	3	0	5	Pass
CT Slope over Z	HU/mm	0.05	0	-0.15	0.15	Pass
CT Standard Deviation over Z	HU	0.3	0.5	0	1	Pass
Contrast to Noise Ratio	HU/HU	2.1	3.2	1	3	Pass

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